# STORZ

## KARL STORZ — ENDOSKOPE

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**1 Information for users of  
KARL STORZ devices**

**It is recommended to check the suitability of  
the products for the intended procedure prior  
to use.**

Thank you for your expression of confidence in  
the KARL STORZ brand name. Like all of our  
other products, this product is the result of years  
of experience and great care in manufacture. You  
and your organization have decided in favor of  
a modern, high-quality piece of equipment from  
KARL STORZ.

This instruction manual forms part of the device.  
KARL STORZ SE & Co. KG, referred to in the  
following simply as KARL STORZ, assumes no  
liability nor provides any warranty whatsoever  
for damage and consequential damages that  
arise due to non-compliance with the instruction  
manual.

- Read the instruction manual carefully before  
using this device.
- Store the instruction manual in a safe place  
throughout the service life of the device.
- Keep the instruction manual accessible to  
surgical staff.
- Give the instruction manual to each successive  
owner and/or user of this device.

**1.1 Revision index device**

| Software version | Last revised      |
|------------------|-------------------|
| Version 1.0.2    | Rev. 03 – 2019/11 |

**1.2 Validity**

This instruction manual Art. no. 96326006D  
applies only to the devices UH 300E and  
UH300UE.

**1.3 Other applicable documents**

- Comply with other applicable documents in  
the appendix.

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## Information for users of KARL STORZ devices

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#### 1.4.1 Structure of warnings

#### **SIGNAL WORD**

Risk type, source and consequences thereof  
(Personal injury)!

- ▶ Measure to avoid the risk

#### **NOTE**

Risk type, source and consequences thereof  
(Material damage)!

- ▶ Measures

### 1.4.2 Risk levels in warnings

| Symbol | Risk level     | Probability<br>of occur-<br>rence | Conse-<br>quences of<br>non-com-<br>pliance |
|--------|----------------|-----------------------------------|---------------------------------------------|
|        | <b>DANGER</b>  | Immediate<br>risk                 | Death,<br>serious<br>injuries               |
|        | <b>WARNING</b> | Possible<br>risk                  | Death,<br>serious<br>injuries               |
|        | <b>CAUTION</b> | Possible<br>risk                  | Minor<br>injuries                           |
|        | <b>NOTE</b>    | Possible<br>risk                  | Material<br>damage                          |

#### 1.4.3 Tips

- Tips and additional information to facilitate tasks

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**1.4.4 Other symbols and labeling**

| Symbol/<br>labeling           | Meaning                                         |
|-------------------------------|-------------------------------------------------|
| ☑                             | Prerequisite for an action                      |
| △                             | Action with one step                            |
| 1.<br>2.<br>3.                | Action with several steps in a binding sequence |
| ☞                             | Result of preceding action                      |
| •                             | List (first level)                              |
| •                             | List (second level)                             |
| <b>Emphasis</b>               | Emphasis                                        |
| ..., see section xxx page xxx | Cross reference                                 |

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### 2 Safety

#### 2.1 Intended use

The high-frequency surgical device is exclusively intended for creating electrical power for high-frequency surgical applied parts.

#### 2.2 Indication(s)

The high-frequency surgical device with the corresponding accessories is suitable for all applied parts which require the energy provided as per the 'Technical data' section of the relevant manual.

#### 2.3 Contraindication(s)

The high-frequency surgical device is not used in direct contact with the patient but rather provides energy and power for HF applied parts. Use with the corresponding applied parts is contraindicated if use of such parts is contraindicated.

The device is not suitable for applied parts which require continuous activation of the HF surgical device.

#### 2.4 Patient profile

The use of AUTOCON® III 300 is not limited to a certain patient profile (sex, age, weight etc.). The patient population is limited by the applied parts connected to the device.

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## 2.5 General safety information

- ▷ Ensure that no electronic devices that are subject to interference from electromagnetic fields are set up in the vicinity of the HF device.
- ▷ Observe the instructions on electromagnetic compatibility provided in section EMC, page 159.

**WARNING:** Danger of electric shock!

- ▷ Always connect the HF device to a mains power system with a protective conductor in order to prevent electric shock.
  - ▷ Ensure constant access to the power plug at the site of installation. The line voltage to the device is only definitively disconnected once the power plug has been unplugged.
- Additional equipment connected to medical electrical equipment must verifiably comply with the respective IEC or ISO standards (e.g. IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see Clause 16 of IEC 60601-1). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements.

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If in doubt, please contact your local specialist dealer or the Technical Support, see section Technical Support, page 68.

- i** To protect personnel, KARL STORZ recommends the use of a smoke evacuator to extract electrosurgical smoke. To this end we recommend the KARL STORZ S-PILOT® UP 501 which, in combination with the AUTOCON® III 300 (plus connecting cable UP 004), guarantees automatic smoke extraction.

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## 2.6 Personal safety instructions

### 2.6.1 Ambient conditions

- Do not use the HF device in the immediate vicinity of the patient. Observe the minimum distances recommended by KARL STORZ (see Fig. 1).

### 2.6.2 Patients with pacemakers

- Malfunction or destruction of the pacemaker can endanger the life of the patient or result in irreversible injuries to the patient.
- In the case of patients with pacemakers, consult the cardiologist before carrying out HF surgery.
- Use bipolar HF methods.
- Attach the HF neutral electrodes close to the operating field.
- Set the demand pacemaker to a fixed frequency.
- Ensure that the pacemaker does not come into contact with the HF electrode.
- Keep a fully operational defibrillator within reach.
- Carry out a postoperative pacemaker check.

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
### 2.6.3 Hazard-free patient positioning – Avoid burns from leakage currents

- ▷ Position patients so that they are not touching any metal parts that are grounded or have considerable capacitance relative to ground (e.g. operating table brackets). If necessary, place anti-static cloths between the patient and the underlay.
- ▷ Ensure that the patient does not touch any wet clothes or underlay.
- ▷ Place anti-static cloths between areas with heavy sweating and skin-to-skin contact areas on the patient's torso.
- ▷ Ensure that the patient is resting on a suitable surface in order to prevent pressure necrosis.
- ▷ Drain urine via a catheter.
- ▷ Do not touch the patient when the surgeon is using the activated HF instrument on the patient.

**i** Due to the defined cable routing and reduced HF voltages and power with bipolar applications there are generally fewer leakage currents than with unipolar applications.  
Yet with some bipolar modes the power and voltages are similar to unipolar modes which can result in higher leakage currents.

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#### 2.6.4 Correct connection of the HF device

- ▷ Always ground the HF device via the potential equalization connector. Also observe the requirements in section 8.6.7 of IEC 60601-1 regarding medical electrical systems.

**i** The device features a potential equalization connector which

must be connected via a potential equalization line with a connector in the room where it is set up. This additional potential equalization measure means that all conductive parts in the vicinity of the patient are connected electrically to each other and small contact voltages are compensated. As such, both transient disturbances from the switching of adjacent devices and the risk to patients from low-frequency electric current are reduced. The connection should be as low-resistance as possible and should not exceed a resistance value of 0.2 Ohm between the central potential equalization of the operating room and the potential equalization connector AUTOCOR® III 300 and every other metallic touchable surface. If the design of the installation does not include an equipotential bonding bar, the connection to the PE conductor should be as low-resistance as possible.

- ▷ With simultaneous use of HF surgical devices and physiological monitoring devices on a patient, the monitoring electrodes should be positioned as far away as possible from the surgical electrodes.
- ▷ Needle electrodes for monitoring are not recommended.
- ▷ In all cases monitoring systems are recommended which include mechanisms for limiting the high-frequency current such as protective resistors and HF chokes.
- ▷ Place lines from monitoring devices so that they do not lie on the patient's skin.
- ▷ Keep the leads to the HF electrodes as short as possible and position them so that they do not touch the patient or other leads.
- ▷ Do not place any objects on the HF device.

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### 2.6.5 Correct use of the HF device

Inadvertent activation of the HF device outside the user's field of vision can injure the patient.

- ▷ Activate the HF device only when the electrode is in your field of vision and you can quickly deactivate the HF device at all times.
- ▷ If the HF device is activated inadvertently, switch it off immediately using the Standby button.
- ▷ Take particular care when using a footswitch or manual switch.

Improper preparation or errors in the HF device can cause damage to the HF device.

- ▷ Use the automatic monitoring functions to ensure that the HF device is working properly. See the section Test for proper functioning, page 29 for information on the automatic test functions.
- ▷ Ensure that no conductive fluids (e.g. blood or amniotic fluid) have penetrated the footswitch or the manual switch.
- ▷ Ensure that the cables for the footswitch and the manual switch are free from short circuits and broken leads.

### 2.6.6 Configuring HF device settings and using accessories

Setting the output power too high can injure the patient! Therefore, before you increase the output power, ensure that:

- the neutral electrode is attached properly;
- the working electrodes are clean;
- the plug connections are all correct.

A clearly lower output value or functional failure of the HF surgical unit during normal operation can be due to insufficient attachment of the neutral electrode or insufficient contact in the connections. In this case, check attachment of the neutral electrode and the connections before selecting a higher output power.

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### Setting the HF device correctly

- ▷ To prevent inadvertent (thermal) tissue damage during operations on body parts with small cross sections and in areas with high resistance (bones or joints), use the bipolar method in these areas.
- ▷ Set the level of the acoustic signal that sounds when the electrode is activated so that it is always clearly audible.

### Nerve and muscle excitation due to low-frequency currents!

During HF surgical operations (especially when an arc is formed), part of the HF current is converted into a low-frequency current. This current can trigger muscle contractions in the patient:

- ▷ To minimize the risk of injury to the patient, set the power and the effect as low as possible.

### Correct use of accessories

- ▷ Only use insulated accessories.
- ▷ Check all electrodes for sharp edges and projecting parts before use.
- ▷ Only use electrodes that are free of defects and in good working order.
- ▷ Never place active electrodes on or near the patient. Ensure that the electrodes cannot come into direct or indirect contact (through electrically conductive objects) with the patient.
- ▷ Electrodes which are not currently in use must not come into contact with the patient.
- ▷ Place the instruments in a safe place: sterile, dry, non-conductive, clearly visible. Instruments which have been placed down must not come into contact with the patient, the medical personnel, or flammable materials.
- ▷ Do not remove hot electrodes from the patient's body directly after cutting or coagulation.
- ▷ Ensure that there is sufficient distance between the patient cables and the cables of the HF device.
- ▷ Do not run the patient cables across the patient.

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## 2.7 Product-related safety instructions

Devices manufactured by KARL STORZ are developed in accordance with the current state of technology and generally accepted safety rules. Despite this, using these products can lead to risks to the life and health of the user or third parties and/or damage to the device or other objects.

- ▷ Only use accessories approved by KARL STORZ, see section Accessories/Replacement parts, page 158.
  - ▷ Only use the device if it is free of technical defects and in good working order and only for the intended purpose, always remaining aware of safety requirements and risks and complying with this instruction manual.
  - ▷ Have faults that can adversely affect safety (e.g. deviations from the permissible operating conditions) rectified without delay.
  - ▷ Wipe down the HF device only with cleaning agents and disinfectants that are approved in the country of use for surface cleaning. Other reprocessing methods are not permitted! (See section Cleaning and disinfection, page 64.)
  - ▷ If any fluids penetrate the device, drain them immediately.
- Damage to the device can lead to an undesirable increase in output power due to a device malfunction.
- Certain devices or accessories can cause danger in lower power settings.

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## 2.8 Safe handling (general instructions)

- ▷ Before each use of the device, check to ensure that it is functioning properly and is in good working order and connected properly.
- ▷ Observe the instructions on intended use in conformance with standards, see section Error indications for EASY monitoring, page 61.
- ▷ During use, always observe and comply with the acoustic signals and/or error indications of the HF device, see section Error indications for EASY monitoring, page 61.
- ▷ The device and accessories may be operated and used only by people who have the necessary training, knowledge and experience.
- ▷ Regularly inspect the accessories, especially the electrode cables, endoscopic accessories and neutral electrodes, for damage to the insulation, proper operation and expiration date.
- ▷ Never place the instruments on the patient or on the devices.
- ▷ Ensure that no instruments are being cleaned when AUTOSTART is activated.
- ▷ Wear suitable gloves during operations.

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### 2.8.1 Operation area: Avoiding ignition and explosions


Sparks are generated when the HF device is used as intended!

- ▷ Do not use the HF device in areas where there is a risk of explosion.
- ▷ Do not use any flammable or explosive liquids.
- ▷ If the display fails, do not use the HF device any longer!
- ▷ During operations (e.g. in the head or thoracic regions), avoid using ignitable anesthetics and gases that support combustion (e.g. nitrous oxide or oxygen) or extract them using a vacuum system.
- ▷ Use exclusively non-flammable cleaning agents, disinfectants and solvents (for adhesives). If you cannot avoid the use of flammable cleaning agents, disinfectants or solvents, ensure that they have fully evaporated before using HF surgical equipment.
- ▷ Ensure that no flammable liquids collect beneath the patient or in body cavities (e.g. the vagina). Suction or flush body cavities before activating the device.
- ▷ Wipe off all liquids before using the HF device.
- ▷ Ensure that no ignitable endogenous gases are present.
- ▷ Ensure that all flammable materials saturated with oxygen (e.g. cotton or gauze) are kept far enough away from the HF environment that they cannot ignite.

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#### Application point for neutral electrode

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#### 2.8.2 Applying the neutral electrode

**i** Observe the instructions on the use of the neutral electrode in the instruction manual and the information on the package of the neutral electrode.

In the unipolar HF method, the neutral electrode feeds the current introduced into the patient's body at the surgical site back to the HF device.

▷ To prevent a rise in temperature at the current exit point, the following conditions must be ensured:

- Sufficiently large contact surface between the neutral electrode and the patient's body;
- High electrical conductivity between the neutral electrode and the patient's body.

▷ To prevent the patient being burned by the neutral electrode, you must comply with the following conditions:

- Select the application point for the neutral electrode so that the current paths between the active and neutral electrodes are as short as possible and run longitudinally or diagonally to the patient's body (because muscles are more conductive in the direction of the fibrils, see Fig. 2).
- For surgery in the thoracic region, do not run the current path transversely across the patient's body and ensure that the patient's heart is never in the current path.
- Depending on the surgical site, apply the neutral electrode to the nearest upper arm or thigh if possible, but never closer than 20 cm.
- In the case of self-adhesive disposable electrodes, comply with any further manufacturer instructions regarding the point of application.
- Ensure that the application area is free of scar tissue, bony protuberances, surface hair and ECG electrodes.

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- Ensure that there are no implants (e.g. bone nails, bone plates or endoprotheses) in the current path.
- Ensure that no short circuits can occur at the neutral electrode connection.
- Avoid areas where fluids may collect.
- As far as possible, use split neutral electrodes with a sufficiently large surface area (patient age and max. output power during operation have to be considered).

#### Before applying the neutral electrode

- ▷ Remove excess body hair.
- ▷ Clean the application site, but do not use any alcohol, since it dries out the skin and increases the transition resistance.
- ▷ If the patient has poor circulation, massage or brush the application site.
- ▷ Attach the neutral electrode over the entire contact surface evenly. Secure reusable neutral electrodes with rubber bands or elastic straps so that they do not loosen when the patient moves. Ensure that the patient's circulation is not impaired (risk of necrosis).
- ▷ Never use wet cloths or conductive pastes.
- ▷ Ensure that no liquids (e.g. cleaning fluids, disinfectants, blood or urine) penetrate between the patient and the neutral electrode.
- ▷ Do not place the neutral electrode under the patient's buttocks or back.
- ▷ Ensure that there are no ECG electrodes in the current path of the HF device.
- ▷ Check the neutral electrode for damage and to ensure that it is working properly.
- ▷ Replace defective accessories immediately.

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### Example application using a disposable electrode

- ▷ Remove the protective film and attach the self-adhesive disposable electrode to the patient. Ensure that the long edge of the disposable electrode faces the operation site and the electrode is fully in contact with the skin. This avoids excessive current concentration on the short edge.
- ▷ Using both hands, press the self-adhesive disposable electrode firmly against the patient's skin.
- ▷ Clamp the electrode tab to the electrode cable.
- ▷ After the operation, remove the disposable electrode carefully to avoid skin damage.

### With a one-piece neutral electrode

- ▷ Check the one-piece neutral electrode during surgery.
- ▷ Ensure that the one-piece electrode is not blocked at the device.

### With a split neutral electrode

- ▷ Apply the split neutral electrode correctly and without any additional objects, as the HF device does not recognize the bridging of the split surfaces with other objects.
- ▷ Ensure that the HF current flows equally to both parts of the split neutral electrode.

**i** For monitoring the neutral electrode, see section EASY neutral electrode monitoring (EASY monitoring), page 31.

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- 17 Connection socket 1 for footswitch
  - 18 Connection socket 2 for footswitch
  - 19 Connection for potential equalization
  - 20 Power supply for IEC connector
  - 21 Fiber optic signal input socket
  - 22 Fiber optic signal output socket
- Use the following connections only for service and training purposes:**
- 23 Ethernet connector (only for connection to KARL STORZ OR1™ control NEO)
  - 24 USB connector
  - 25 Connector S-PILOT® UP 501 with UP 004 S-PILOT® connecting cable, dia. 3.5 mm, 300 cm (max. 25V/1 A)
  - 26 UART communication interface
  - 27 Power switch
  - 28 Line fuses (see page 70)

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-  The Ethernet interface is only suitable for connection of the OR1™ control NEO for a dedicated network.
-  Only use the USB interface as described in this manual. To this end, only use the KARL STORZ USB stick 20 0402 82.

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### 3.2 Symbols used on the device

| Symbol                                                                                | Designation                                                                                                                                        |
|---------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Footswitch connection                                                                                                                              |
|    | Neutral electrode isolated from ground for HF                                                                                                      |
|    | CF type applied part with defibrillation protection (cardiac floating)                                                                             |
|    | Alternating current                                                                                                                                |
|    | Standby                                                                                                                                            |
|   | During activation (of the HF device) HF energy in the radio frequency range 9 kHz to 400 GHz is applied, which produces electromagnetic radiation. |
|  | This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE), see Disposal.           |
|  | Labeling of the (active) HF output; Caution: Hazardous voltage.                                                                                    |

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| Symbol                                                                                | Designation                                                                                                                     |
|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|
|    | Manufacturer                                                                                                                    |
|    | Date of manufacture                                                                                                             |
|    | Follow instructions for use                                                                                                     |
|    | Potential equalization connector                                                                                                |
|    | Fiber optic signal input                                                                                                        |
|    | Fiber optic signal output                                                                                                       |
|    | Ethernet connector<br>(connection to KARL STORZ OR1™ control NEO system).                                                       |
|    | USB connector*                                                                                                                  |
|    | S-PILOT® connector*<br>(surgical smoke extraction)                                                                              |
|    | UART communication interface for service purposes*                                                                              |
|    | Reference number                                                                                                                |
|  | Serial number                                                                                                                   |
|  | The device satisfies the Medical Device Directive (MDD) 93/42/EEC (code number 0123 indicates notified body).                   |
|  | Consult instructions for use                                                                                                    |
|  | General warning of physiological effect: The device generates physiological effects which may endanger the patient and/or user. |
|  | 'Rx only'                                                                                                                       |

\* In the event of queries contact Customer Service (see section 2.3).

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### 3.2.1 Symbols on the packaging

For the meanings of the symbols printed on the label or packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no.96216316DF. It can be downloaded from [www.karlstorz.com](http://www.karlstorz.com).

### 3.3 Assignment of color to function (on device and accessories)

| Color  | Function    |
|--------|-------------|
| Yellow | Cutting     |
| Blue   | Coagulation |

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### 3.4 Scope of supply\*

1 AUTOCON® III 300 UH 300E  
1 Power cord 400 A

1 Instruction manual with inspection protocol

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### 3.5 Components required for operation

- Power cord
- Footswitch
- Neutral electrode for unipolar applications
- Connecting cable for neutral electrode or instrument
- Instrument (unipolar or bipolar)

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### 3.6 Operating conditions

Temperature: +10°C to +40°C  
Relative humidity: 30 to 75%, non-condensing  
Air pressure: 700 to 1 060 hPa  
Operating altitude (max.): 3000 m above sea level

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### 4 Preparation

#### 4.1 Setting up the HF device

**NOTE:** HF equipment generates high-frequency electromagnetic fields that can interfere with the operation of sensitive electronic equipment and, for example, lead to image distortion. HF equipment, and particularly its cabling, should thus never be installed close to any equipment sensitive to electromagnetic interference. Parallel cable routing with cables from other devices should be avoided as a general rule, as HF cables emit HF emissions when activated, which may interfere with other equipment. Additional measures include:

- Connect the power connection of the AUTOCON® III 300 or the impaired device to another circuit.
- Increase the distance between the AUTOCON® III 300 and the impaired device.
- Lay instrument cables so that, as far as possible, they are not in the vicinity of other devices and their connecting cables.

**WARNING:** Risk of electric shock! Always connect the HF device to a mains power system with a protective conductor in order to prevent electric shock.

**DANGER:** Risk of burns to patients due to excessive leakage current! Locate the HF device outside the immediate vicinity of the patient, see section Ambient conditions, page 6.

HF devices may be used only in rooms used for medical purposes that meet the requirements of DIN VDE 0100-710 or IEC 60364-7-710.

If the HF device was previously stored or transported at temperatures below +10°C or a relative humidity above 75%, non-condensing, it will take approximately three hours to adjust to room temperature.

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#### 4.2 Switching on the HF device

**i** Do not use the HF device if the display components are not working! For troubleshooting see section Recognizing and remedying errors, page 53.

1. Switch the HF device on using the power switch on the rear side of the device and then touch the permanently lit Standby button on the front panel.

The HF device performs a self-test: All display and control elements light up.

2. Check all display and control elements for proper operation:
  - Power switch
  - Touch screen
  - Unipolar connection sockets
  - Bipolar connection sockets
  - Activation bar for unipolar and bipolar sockets

The main screen appears, and the HF device is ready for use.

The parameters of the most recently selected program appear on the display.

**i** If the device is switched off at the front panel, the Standby mode is activated.

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#### 4.3.1 Instruments for unipolar application

1. Plug the neutral electrode cable into the socket for the neutral electrode and choose the corresponding neutral electrode type, see chapter Selecting the neutral electrode, page 39.
  - ☛ The neutral electrode button will change from gray to the measurement colour (green, yellow or red).
2. Connect the electrode handle to the unipolar connection socket (see Fig. 3).
  - or –
 In the case of an accessory without a finger button: Connect the footswitch and unipolar connecting cable.
  - or –
 Connect the unipolar cable to the unipolar connection socket for unipolar instruments.

#### Unipolar connection socket ⑨ (see Fig. 3)

- Connection socket ⑨ 'top':
  - ⑪ 3-pin US type
  - ⑫ 5-mm KARL STORZ/Erbe VIO
  - ⑬ 4-mm socket (foot-switch)

#### Connection socket for neutral electrode (see Fig. 4)

- ⑩ The peak output voltages of the different modes are listed in chapter 10.2 (pages 97-134).

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#### 4.3.2 Instruments for bipolar application

1. Connect the bipolar cable to the instrument (e.g. forceps).
  2. Connect the bipolar cable to the bipolar connection socket (see Fig. 5).
  3. For bipolar use without AUTOSTART, connect the footswitch.  
– or –  
Select a mode with AUTOSTART function at the corresponding socket.
- ☛ Once the instrument is connected, the application starts after the configured reaction time.

**i** When introducing bipolar instruments into a metal trocar, the AUTOSTART function should be disabled beforehand!

#### Bipolar connection sockets

- Bipolar socket ⑭ (see Fig. 5)  
⑮ 2-pin US type (28.58 mm)  
⑯ 2x KARL STORZ/Erbe VIO
- \* Applied part type CF according to IEC 60601-1, with defibrillation protection

**i** The peak output voltages of the different modes are listed in chapter 10.2 (pages 135-156).

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
### 4.3.3 Connecting the footswitch

In addition to the manual switch, a footswitch can be used to activate various operating modes.

▷ Connect the desired footswitch only during operation to one of the two connection sockets for footswitches.

☛ The HF device automatically detects the connected footswitch and indicates this on the front panel display, including the selected connection socket.

❗ One dual-pedal footswitch and one single-pedal footswitch can be connected. Footswitches without a changeover switch cannot be used (see table).

During the operation, in addition to the footswitch, only, if relevant, the connection to the KARL STORZ OR1™ control NEO system may be plugged into the rear of the device.

The following footswitches can be connected to the HF device:

| Art. no. | Designation                                             |
|----------|---------------------------------------------------------|
| UF 901   | Single-pedal footswitch with changeover function switch |
| UF 902   | Dual-pedal footswitch with changeover function switch   |

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#### 4.4 Test for proper functioning

##### 4.4.1 Automatic test function

The HF device automatically performs cyclic testing during operation. If any errors occur, see section Recognizing and remedying errors, page 53.

##### 4.4.2 Performing the test for proper functioning

Perform the following functional test before putting the device into service:

- ❗ The accessories must be designed for the specified maximum voltage.

1. Connect the neutral electrode and attach it securely to the patient's arm.

- ➡ The EASY neutral electrode indicator changes to green.

2. Remove the neutral electrode.

- ➡ The indicator changes to red, acoustic signals sound.

- ❗ The neutral electrode used for this test may not later be used for an operation.

3. Connect a unipolar HF handpiece to a unipolar socket if there is a green EASY indicator and use the manual switch and footswitch to individually activate 'Cut' and 'Coag'.

4. Check the settings on the display.

5. Now change to the bipolar output and connect bipolar forceps.

6. Select a mode with AUTOSTART, grasp moist gauze with the forceps, and check the display.

7. Now change to a mode without AUTOSTART and use the footswitch to activate the bipolar output. Check the settings and indicators in the bipolar section.

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#### 4.4.3 Actions in case of problems

Proceed as follows in case of functional problems:

1. Immediately disconnect the patient from the HF device.
2. Perform a technical inspection of the HF device.
3. Observe the relevant national regulations. For example, in Germany report incidents and near-accidents to the German Federal Institute for Drugs and Medical Devices in accordance with Section 3 of the German Ordinance on the Installation, Operation and Use of Medical Devices (MPBetreibV).  
Observe the provisions of the in-house reporting system in this regard.
4. Consult the Technical Support, see section Technical Support, page 68.

The HF device can be switched off at any time using the power switch ② on the rear of the device.

#### 4.5 Neutral electrode monitoring

Always use the largest possible electrode when attaching a neutral electrode.

##### 4.5.1 General information

KARL STORZ recommends using split neutral electrodes, since only this type of electrode allows the HF device to detect detachment of the neutral electrode if this occurs.

Monitoring of the neutral electrode minimizes the risk of burns at the site where the neutral electrode is attached.

Two types of neutral electrodes can be monitored:

- Split small\* neutral electrodes (for use with reduced power)  
\* available commercially as baby or children electrodes
- Split neutral electrodes

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The type of neutral electrode and its contact quality are selected and/or shown in the Neutral electrodes menu. Errors related to the neutral electrode and possible remedies are shown on the display, see section Recognizing and remedying errors, page 53 (see Fig. 7, neutral electrode modes).

#### 4.5.2 EASY neutral electrode monitoring (EASY monitoring)

The maximum output power of unipolar current types is reduced to 50 W when a small ('children') electrode is selected.

The EASY monitoring function measures changes in the resistance between the patient and the high frequency surgical unit before and during HF activation. If necessary, it generates visual and audible alarms to request staff intervention. This requires using a split neutral electrode with appropriate contact areas and suitable transition resistance, attached to the patient according to the manufacturer's instructions. The EASY system does not monitor the currents through the individual contact surfaces of the split neutral electrode.

A neutral electrode with a surface area of at least 90 cm<sup>2</sup> must be used for the unipolar modes 'Resection' and 'Contact Soft'.

If an error message is generated, the display changes from green via yellow to red, depending on the type of error.

**NOTE:** Incorrect application of the neutral electrode! Ensure compliance with the specifications for correct attachment of the neutral electrode with regard to size, adhesive properties and full-surface contact of the complete electrode.

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
## 5 Operation

### 5.1 Connecting the power supply

- ✓ The line voltage must match the voltage specified on the type plate.
- ✓ The potential equalization cable must be connected to the potential equalization connector (see section 4.1). The device's ground line should be installed by a qualified electrician.
- 1. Connect the power cord to the generator and the plug to the grounded wall socket.
- 2. Use the power switch 27 on the rear to switch the device on (see Fig. 8).
- 3. Switch the device on at the Standby button.
- ☛ The device performs a self-test.
- ☛ The bars of the activation indicators light up (see page 17).
- ☛ Full functionality of the loudspeaker is indicated by the start melody.
- ☛ The individually configurable start screen appears if it has been set up.

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## 5.2 Program overview

### 5.2.1 Display

The status bar is located at the top of the display (see Fig. 9). The settings of the two sockets are shown below the status bar. Settings can be configured for each of the sockets.

1. The 'Effect' button is used to set the effect of electro-surgical cutting or coagulation.
2. Assign the footswitch using the 'Pedal' button.
3. The 'Mode' button allows the desired type of current to be selected.
4. Tap on 'max. Watt' to adjust the maximum output power.

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

### 5.3 Activating and deactivating sockets

- ▷ To activate a deactivated socket, plug a connecting cable into the connection socket.
- or –

Press the button 'activate socket' (see Fig. 11).

- ☛ An overview of the socket settings appears.

The overview is grayed out if no instrument is connected to the socket (see Fig. 12).

The socket illumination extinguishes and the socket selection window lights up when an instrument is plugged in.

- ▷ To hide sockets, tap on the 'off' button next to the socket setting overview (see Fig. 13).

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#### 5.4 Unlocking the screen

The device screen locks automatically. To unlock it, tap on any control element and push the slider which then appears from left to right. An open lock then appears in the status bar (see Fig. 14, Unlocking the screen).

To deactivate the key lock or to change the time up to automatic screen lock, see section Menu 'System settings', page 44.

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
15

## 5.5 Configuring output currents

- i** All selection windows are closed after 30 seconds without assumption of the modification.
- i** If selection windows are open, all buttons outside this field are inactive. In this case, activations are not possible.
- i** A change to the currently loaded program, e.g. by adjusting the output, is indicated by the information 'changed' below the program name.

### 5.5.1 Selecting the mode


1. To select the type of current for unipolar cutting, select the socket on the left side (see Fig. 15).
2. Press the yellow icon under the 'Mode' button.  
**i** A selection screen appears for the available modes.
3. Select the desired mode by pressing the corresponding button.  
- or -  
Deactivate the mode by pressing the button 'Cutting off (X)'.
4. For additional selection options press ◀/▶ inside the selection window.
5. Confirm your selection by pressing the 'OK' button.  
**i** The main screen will be displayed.  
- or -  
Press the 'Back' button to return to the main screen without changing the selection.
- i** If a mode is changed within a socket, the set parameters, e.g. effect and max. Watt, remain the same for the respective mode. However, if, for example, the factory default setting of a mode is adapted and subsequently changed to a different mode and then back again, the user changes are not undone.

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
#### 5.5.2 Specifying power limits

1. To select the maximum output for unipolar cutting, first select the unipolar socket and then press the yellow icon under the 'max. Watt' button (see Fig. 16).
2. Use the '+' and '-' buttons to adjust the power level in single steps.  
- or -  
Use the slider to set the power level in steps of 10.
3. Press the '?' button for more information on this selection.
4. Confirm your selection by pressing the 'OK' button.  
- or -  
Press the 'Back' button to return to the main screen without changing the selection.

#### 5.5.3 Selecting the effect

1. To select the effect for unipolar cutting, first select the unipolar socket and then press the yellow icon under the 'Effect' button (see Fig. 17).
2. Use the '+' and '-' buttons to adjust the effect in individual steps.  
- or -  
Use the slider to set the effect.
3. Press the '?' button for more information on this selection.
4. Confirm your selection by pressing the 'OK' button.  
- or -  
Press the 'Back' button to return to the main screen without changing the selection.

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19

### 5.5.4 Assigning the footswitch

Handpieces and instruments with manual switches can be activated without a configuration setting.

A dual-pedal footswitch and/or single-pedal footswitch, each with a changeover switch, can be connected.

The changeover switch enables switching between pedal levels.

1. Call up the footswitch menu by pressing the button 'Pedal' (see Fig. 18).

2. Select the desired footswitch by pressing the corresponding button. For example, choose the active pedal level for cutting and coagulation for the upper left socket.

- or -

Deactivate the footswitch by pressing the button marked with an 'x'.

The edge of the selected button lights up green (see arrow in Fig. 18).

3. Confirm the selection by pressing the 'OK' button.

- or -

Press the 'Back' button to return to the main screen without changing the selection.

The socket is assigned to the active pedal level. Pedal levels can be changed using the footswitches. Press the black changeover switch to change the socket.

The orange background of the changeover switch in the 'Pedal' button indicates that the lower left-hand socket is activated (see Fig. 19).

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**i** The maximum output power of unipolar current types is reduced to 50 W when a small\* electrode is selected.  
\* available commercially as baby or children electrodes

2. Select the type of connected neutral electrode by pressing the corresponding symbol.

EASY: for monitoring split neutral electrodes (see Fig. 21)  
BABY: for monitoring split small neutral electrodes (see Fig. 22)  
MONO: to select a one-piece neutral electrode (see Fig. 23)

3. Press the '?' button for more information on this selection.
4. Confirm your selection by pressing the 'OK' button.

- or -  
Press the 'Back' button to return to the main screen without changing the selection.

**i** The selected type of neutral electrode in connection with a color indicator for the contact quality is shown in the status bar.

**i** When using the 'EASY' and 'BABY' mode, no unsplit electrodes are accepted. Using the 'MONO' mode, no split electrodes are accepted.

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According to the contact quality, the icons are shown in color for neutral electrodes:

| Icon/Button                                                                          | Description                                                                    |
|--------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|
|    | Split neutral electrode contact<br>quality OK                                  |
|    | Split neutral electrode contact<br>quality not optimum                         |
|    | Split neutral electrode contact<br>quality insufficient                        |
|    | Split neutral electrode not<br>connected                                       |
|    | Split small 'Baby' neutral<br>electrode<br>contact quality OK                  |
|    | Split small 'Baby' neutral<br>electrode<br>contact quality not optimum         |
|    | Split small 'Baby' neutral<br>electrode<br>contact quality insufficient        |
|   | Split small 'Baby' neutral<br>electrode<br>not connected                       |
|  | Non split neutral electrode<br>contact quality OK                              |
|  | Non split neutral electrode not<br>detected or contact quality<br>insufficient |
|  | Non split neutral electrode not<br>connected                                   |
|  | Display contact quality                                                        |

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## 5.6 Mode overview

An overview of the current types that can be executed with the HF device is shown below.

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### 5.6.1 Unipolar modes

| Cutting mode symbol                                                                  | Designation         |
|--------------------------------------------------------------------------------------|---------------------|
|    | Cut reduced 200W    |
|    | Cut                 |
|    | Cut reduced 50W     |
|    | Resection           |
|    | Cut 300-400W        |
|    | Cut blend           |
|   | Intermittent 400W 1 |
|  | Intermittent 400W 2 |
|  | Intermittent 400W 3 |
|  | Intermittent 300W 1 |
|  | Intermittent 300W 2 |
|  | Intermittent 300W 3 |

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| Coagulation mode<br>symbol                                                         | Designation          |
|------------------------------------------------------------------------------------|----------------------|
|  | Forciert Coag 1.8kVp |
|  | Contact Soft         |
|  | Forced Coag          |
|  | Resection            |
|  | Spray                |
|  | Forced mixed         |
|  | Forced Coag with Cut |

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**5.6.2 Bipolar modes**

| Cutting mode<br>symbol                                                               | Designation      |
|--------------------------------------------------------------------------------------|------------------|
|  | Bip. cutting     |
|  | Bipolar scissors |
|  | VAP CUT          |

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| Coagulation mode symbol                                                             | Designation      |
|-------------------------------------------------------------------------------------|------------------|
|   | RoBl®            |
|   | Standard plus    |
|   | Standard         |
|   | VAP COAG         |
|   | Standard AUTO    |
|   | Bipolar scissors |
|   | Micro            |
|  | Forced           |

The information and data regarding settings, application points, application duration and instrument use are based on clinical practice. However, these are only basic guidelines which must be tested for suitability by the operator. Depending on individual conditions, it may be necessary to deviate from the provided data. Medical practice is continuously evolving as a result of R&D and clinical experience. This may also make deviations from the provided data necessary.

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

**5.7 Menü**  
Basic settings can be changed, e.g., operational language, audio, display and saving options.

**5.7.1 Overview**  
The functions shown in the figure on the left (Fig. 24) are available.

**Selecting a menu**  
Press the corresponding selection button to open the menu.

**Exiting a menu**  
Press the 'back' button to return to the main screen.

**5.7.2 Menu 'System settings'**  
In the 'System settings' menu the settings shown in the figure on the left (Fig. 25) can be changed.

**Autostart delay**  
The time can be set under Autostart delay after which the HF generator is automatically activated impedance-dependent following tissue contact (see Fig. 26).

Incremental setting  
Press the buttons '+' and '-'.  
Fast setting  
Move the slider to and fro.

**Language selection**  
The following languages can be selected under 'Language selection':  
German, English, Spanish, French, Portuguese, Italian, Russian, Polish, Turkish, Japanese, Korean, Thai, Indonesian, Chinese, Czech, Arabic, Hungarian, Swedish

**Key lock**  
In order to switch off the automatic screen lock or to set its duration press 'key lock'. The duration can be set between 30 seconds and 5 minutes.

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### OR system selection

The connection to the KARL STORZ OR system OR1™ can be created under the menu item 'OR system selection' (see Fig. 27). Upon delivery the KARL STORZ OR system OR1™ is in Standby mode. This is evident by the symbol .

To create a connection to your OR1™ a LAN cable specified by KARL STORZ must be used. You can get more detailed information on this from our service team. If the OR system OR1™ is in Standby mode, a connection is established automatically when the connecting cable is inserted.

### Deactivating the OR system

(operating room without OR1 system):

1. Tap on 'off'.
2. Press OK to confirm.
- ☛ The OR1™ symbol is hidden. An automatic connection upon connecting the OR1™ system is no longer possible.

Returning the OR system (after deactivation) to Standby mode

1. Select 'OR1, Karl Storz'.
2. Press OK to confirm.
- ☛ The OR1™ symbol is displayed .

### OR1 configuration

(only for trained service personnel)

Special network configurations are possible by tapping on the symbol  (see Fig. 27).

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### 5.7.3 Menu 'Volume'

Use the 'Volume' menu to set the volume of the individual acoustic signals (see Fig. 28).

#### Incremental setting

- ▷ Press the buttons '-' and '+'.

#### Fast setting

- ▷ Move the slider to and fro.

The volume of the activation signal should be increased as necessary for use in relatively noisy surroundings. The alarm tones have a minimum volume and limited changeability.

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| Mode                  | Category         | Frequency (Hz) | Signal type      |
|-----------------------|------------------|----------------|------------------|
| Unipolar Cut          | Activation tones | 635            | Continuous sound |
| Unipolar Coag         | Activation tones | 475            | Continuous sound |
| Bipolar Cut           | Activation tones | 565            | Continuous sound |
| Bipolar Coag          | Activation tones | 505            | Continuous sound |
| Footswitch changeover | Activation tones | -              | Signal tone      |
| ZAP mode              | Activation tones | -              | Signal tone      |
| Fault                 | Alarm tones      | -              | Signal tone      |
| Warning               | Alarm tones      | -              | Signal tone      |
| Note                  | Alarm tones      | -              | Signal tone      |

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### 5.7.4 Menü »Service«

Service contact data can be called up in the 'Service' menu (see Fig. 29).

After entering a password, you can make further settings, e.g. resetting to the delivery state. To access the service level ('Service Tools') (see Fig. 30), enter the password 001224.

### Saving device settings

With the function '1 Backup device' device settings can be saved on the KARL STORZ USB stick 20 0402 82. This includes all saved programs and system settings.

### Installing device settings

Use '2 Restore device' to transfer saved device settings on the KARL STORZ USB stick to an AUTOCON® III 300.

### Adding a start screen

Select the function '3 Add logo' to display a user-defined logo when starting up.

### Deleting a start screen

Select the function '4 Remove logo' to remove a user-defined logo.

### Saving logfiles

Select the function '6 Save logfiles' to save all logfiles onto the KARL STORZ USB stick.

### Resetting to default setting

Select the function '7 Reset to default' to reset all settings and programs to the default setting.

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### 5.7.5 Menü »System information«

In the 'System information' menu various system parameters such as version, serial number, TSI date for AUTOCON® III 300 as well as options (see Fig. 31) are displayed.

For more details see the section Safety inspection, page 65.

### 5.7.6 Menu 'Programs'

1. In the 'Programs' menu you can select programs from a list and assign to or delete from favorites.

- or -

Tap on the current program name in the main screen to access this menu quickly.

2. Press the relevant program name to select the program.
3. To navigate vertically in the program list, press the corresponding button  $\blacktriangle$  /  $\blacktriangledown$ . The programs are always sorted in alphabetical order.
4. Assignment to favorites is possible using the 'star buttons' at the bottom of the screen. Add programs to favorites using the 'green arrow' star button  $\star$  (see Fig. 32 or 33 – left column) and remove these using the 'red arrow' star button  $\star$  (see Fig. 34). The button  $\star$  appears in place of the  $\star$  button if a program is marked in the left-hand column.
5. Sort your favorites using the  $\star$  /  $\star$  star buttons.
6. Load the selected program by pressing 'OK'.

- or -

Return to the main screen by pressing 'back'.

- A maximum of 200 favorites can be selected.

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Saved programs can be deleted in the 'Programs' menu.

1. To this end, select the program to be deleted in the program list by tapping on the program name (see Fig. 32 or 33). To navigate vertically in the program list, press the corresponding button /.
2. To delete the selected program permanently, select the 'recycle bin' button.  
 The selected program is deleted upon confirmation of a prompt.

### 5.7.7 Menu 'Favorites'

The defined favorites can be selected in the 'Favorites' menu. A fast selection of the favorites is possible using the star button in the main screen. Access the next page of the favorites list using the buttons ◀/▶ in the lower part of the screen.

1. Confirm with 'OK' to accept the selection.
2. Press the 'back' button to return to the main screen.

### 5.7.8 Menu 'Save program'

In the 'Save program' menu (see overview under 5.7.1) the current setting can be saved under the same or another program name (see Fig. 35). This function is enabled as standard (upon delivery).

Program names can be created using the keyboard. Symbols, capital or small letters or numbers are selection options.

Program names spread over two lines can be assigned using the 'Enter' button.

1. Confirm with 'OK' to accept the selection.
2. Press the 'back' button to return to the main screen.

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### 5.7.9 Socket settings

In the 'Socket configuration' menu the settings shown in the left-hand figure (Fig. 36) can be changed.

#### Expanding sockets

1. Call up 'Expand socket' in the 'Socket configuration' menu (see Fig. 37).
2. Tap on  next to the two connectors to change the current setting (divide the bipolar socket of the device into 2; see page 18, number ④).

With active socket expansion the  on the left of the two connectors changes to .

3. Confirm with 'OK' to accept the selection.
  4. Press the 'back' button to return to the main screen.
- You can now use the bipolar socket on the bottom of the device as two sockets.

### ZAP mode

With the ZAP mode you can switch between two preset current forms for the same instrument.

1. Call up the ZAP mode in the 'Socket configuration' menu (see Fig. 38).
2. Activate or deactivate the changeover for the individual sockets by pressing the 'tick' next to the relevant socket (see Fig. 39).

If the ZAP mode is active, the 'tick' is displayed filled in.

Confirm with 'OK' to accept the selection.

Press the 'back' button to return to the main screen.

In the main screen a second level of the socket settings is shown for the active ZAP mode (see Fig. 40).

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
You can switch in the main screen using either the handpiece or foot pedal.

For manual switching in the main screen tap on the ZAP mode symbol next to the socket settings (see Fig. 39).

On changing the ZAP mode level, the color of the symbol changes from black to white (see Fig. 40).

Another way of changing the ZAP mode level is by actuating both buttons on the handpiece.

You can also change levels via the black changeover switch on the footswitch.

▷ To do this, select the ZAP mode under the 'Pedal' icon (see Fig. 41).

The ZAP mode symbol now appears in the main screen below the 'Pedal' icon (see Fig. 42).

▷ To change the ZAP mode levels, actuate the black switch on the footswitch.

### 5.7.10 Menu 'System messages'

The 'System messages' menu can be accessed in the main screen via the 'Information' symbol.

Here the errors saved in the HF device which have occurred since switching on the HF device are called up (see Fig. 43).

These messages are not saved when switching off the HF device.

#### Calling up saved system messages:

1. Select a system message.
2. The selected system message is displayed again with '?'.
3. Press 'OK' to return to the overview.

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### 5.7.11 Procedures

Procedures are programs in which parameters such as voltage and output can be determined for specialist medical interventions. You can create your own procedures/programs (see section 5.7.8).

- i** The procedure name is displayed in the status bar (see section 5.2.2).

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
### 5.7.12 Switching off the HF device (decommissioning)

1. Switch off the device at the Standby button on the front panel.
2. Switch off the device at the power switch **27** on the rear (see Fig. 44).
3. Remove connected accessories from the device and, if necessary, reprocess.

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- System error
- EASY monitoring error

### 6.1 System information

System information is displayed as a message on the screen (see Fig. 45).

For further information on the cause and how to remedy this message, press the 'i' button in the status bar.

System information can be divided into three categories:

- Notices (gray display)
- Warnings (orange display)
- Faults (red display)

Notices are shown for five seconds on the screen. Warnings and faults are shown for ten seconds.

While a fault is present, activations are prohibited.

This information can be called up again under the now orange-framed 'i' in the status bar.

These messages are also available in the menu system information and will be deleted when switching off the device.

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The following table describes the cause of the error and the appropriate corrective action.

| Header                         | Message text                                                                                                                                                                                                          |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Confirmation of AUTOSTART mode | You have selected an AUTOSTART mode.<br>Setting the AUTOSTART mode can result in unintentional coagulations, e.g. when bipolar forceps are used for gripping while the AUTOSTART mode is on.                          |
| AUTOSTART Fault                | The instrument is in contact with tissue.<br>AUTOSTART cannot be selected when the instrument is in contact with tissue. Open up the instrument.                                                                      |
| TSI Warning                    | The annual Technical Safety Inspection (TSI) is due.                                                                                                                                                                  |
| Neutral Electrode Fault        | No neutral electrode connected. No neutral electrode is connected. Connect a neutral electrode.                                                                                                                       |
| Neutral Electrode Fault        | Wrong neutral electrode connected.<br>The selection does not match the neutral electrode connected. Connect the neutral electrode which matches the selected mode, or change the mode to match the neutral electrode. |
| Neutral Electrode Fault        | Poor contact with the patient.<br>The resistance between the neutral electrode and the patient is too high. Establish better contact of the neutral electrode.                                                        |
| Neutral Electrode Warning      | Poor contact with the patient.<br>The contact resistance between the neutral electrode and the patient is getting worse. Establish better contact of the neutral electrode.                                           |
| Neutral Electrode Fault        | No cable for neutral electrode connected.<br>No neutral electrode cable is connected. Connect a neutral electrode.                                                                                                    |
| Neutral Electrode Notice       | No cable for neutral electrode connected.<br>The cable for neutral electrode has been removed. Unipolar activation is not possible.                                                                                   |
| Mode Fault                     | No mode selected. No mode was selected for this type of activation. Select the desired mode or change the footswitch assignment.                                                                                      |
| Mode Fault                     | This mode is not allowed for baby neutral electrodes. Use split neutral electrodes with a large surface and the suitable mode for this.                                                                               |
| Mode Fault                     | This mode is not allowed for this socket. The current mode remains active.<br>Choose another socket for this mode.                                                                                                    |
| Footswitch Fault               | No compatible footswitch. The connected footswitch is not compatible with this device. Connect a compatible footswitch with a changeover switch.                                                                      |
| Footswitch Fault               | Footswitch not assigned to a socket.<br>The footswitch has not been assigned to a socket yet.<br>Assign a socket to the footswitch using the 'Pedal' button.                                                          |
| Footswitch Fault               | Error on footswitch connection.<br>Check the footswitch. If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.                                                                |

| Header                                    | Message text                                                                                                                                                                                                                                                                                                         |
|-------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Fingerswitch Fault                        | Fault on fingerswitch connection.<br>Check the handpiece and the connection cable. Please replace them if damaged. If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.                                                                                                     |
| Temperature Warning                       | The temperature of the device is higher than normal.<br>The temperature of the device is elevated. This leads to a reduction of the maximum power.                                                                                                                                                                   |
| Limitation of Continuous Activation       | The maximum activation time has been exceeded.<br>Please only activate the generator in short intervals, in order to avoid harming the patient and damaging the connected instruments or the generator.                                                                                                              |
| Activation Fault                          | While switching on the device, there is an activation by footswitch, fingerswitch or AUTOSTART.<br>Check the handpieces or footswitches for malfunctions.<br>Disconnect the handpieces / footswitches from the device. If the error persists, please contact the Technical Support. Contact: MENU – SERVICE.         |
| Activation Fault                          | There is an activation while connecting the footswitch or fingerswitch.<br>Check the handpieces or footswitches for malfunctions.<br>Disconnect the handpieces / footswitches from the device. If the error persists, please contact the Technical Support. Contact: MENU – SERVICE.                                 |
| Activation Fault                          | There is no instrument connected on the activated socket.<br>Connect an instrument on the designated socket.                                                                                                                                                                                                         |
| Activation Warning                        | The mode for safety inspections is active. Activation is not possible. Quit this mode before activating again.                                                                                                                                                                                                       |
| Intermittent 300W/400W Warning            | Polypectomy snare not in contact with tissue, or check connection cable at snare or generator.<br>Please apply the snare and reactivate. First of all establish contact between tissue and polypectomy snare, or check the connection cable at the snare or the generator. Then activate with the yellow foot pedal. |
| Internal Error xxxx (e.g. with xxxx=4183) | If this message appears again, please contact the Technical Support. Contact: MENU – SERVICE.                                                                                                                                                                                                                        |

With internal errors a number is shown in the header.  
Please quote this number to the Technical Support.

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## 6.2 Error indications for EASY monitoring

Error indications are displayed in three stages (green, yellow and red) when problems occur. When working with a split neutral electrode, the following errors may occur:

| EASY monitoring                       | Cause                                                                                                               | Indication                                                             | Remedy                                                                                                                                                                           |
|---------------------------------------|---------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Lights up yellow                      | Significant increase in resistance<br>Depending on the indication, there may be heating under the neutral electrode | –                                                                      | Stopping the application is not necessary.<br>Check the positioning of the neutral electrode.                                                                                    |
| Switches from green to continuous red | A significant problem occurs when the unipolar current is activated.                                                | An acoustic signal sounds.<br>An error message appears on the display. | Check the neutral electrode and neutral electrode cable, see section EASY monitoring, page 31.<br>▷ Check the neutral electrode cable for proper connection and external damage. |
|                                       | Loosened electrode                                                                                                  | An acoustic signal sounds.<br>A warning message appears on the display | ▷ Correct the positioning of the neutral electrode. If the error persists, replace the neutral electrode.                                                                        |

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## 7 Reprocessing

**!** **WARNING:** Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

**!** **WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

**!** **CAUTION:** When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

**!** **CAUTION:** National laws and regulations must be observed.

**i** The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting [www.karlstorz.com](http://www.karlstorz.com).

### 7.1 Reprocessing of accessories

▷ The reprocessing of accessories is described in the respective manual.

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### 7.2 Cleaning and disinfection

**NOTE:** Damage to the HF device due to incorrect reprocessing: No method other than the wipe-down disinfection described here is permitted.

**WARNING:** Fire hazard and risk from electric shock! Disconnect the plug before cleaning. Only use approved cleaning/disinfection agents for surface cleaning as per the manufacturer's instructions. It is essential to ensure that no liquids can penetrate the housing. Make sure that the AUTOSTART function is deactivated.

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### Wipe-down disinfection of devices and footswitches

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry and low lint cloth.

### Reprocessing limits

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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## 8 Maintenance and repair

### 8.1 Maintenance

**DANGER:** Risk of infection! Carry out a surface disinfection and wrap the device in addition to the shipping packaging material before allowing the device to leave the hospital or practice to avoid spreading germs and infections.

- ▷ Check the device, the equipment cart and the accessories (e.g. footswitch, cable) after each use for damage or defects. In particular, make sure that the insulation is intact on all cables.
- ▷ Do not use any damaged device, damaged equipment cart or damaged accessories.
- ▷ Replace defective accessories immediately.
- ▷ Have the safety inspection for the device performed once a year. Please consult and comply with the respective service manual for additional technical information.

### 8.1.1 Safety inspection

Safety inspections must be performed once a year.

- ▷ The next safety inspection date of the AUTOCON® III 300 can be displayed in the menu, see section Menu 'System settings', page 48.

A warning message appears during system start-up if a safety inspection is due. Press OK to confirm this message.

Any shorter safety inspection cycles specified in national regulations must be observed.

- ▷ The device and accessories may be inspected only by persons who have the required training, knowledge or experience and who can perform the inspection independently.

- ▷ With regard to the safety inspection, you must comply with the country-specific rules and regulations.

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The inspector documents the inspection results and measured values corresponding to the printed inspection protocol.

In the case of severe deviations from the values of the attached final acceptance report, or if the specified maximum values were exceeded:

- ▷ Send the HF device to the service center, see section 9.1 – Technical Support, page 68.

## 8.2 Repair

**NOTE:** You can damage the HF device by doing your own repairs and modifications of medical equipment!  
If a repair is necessary, only allow this to be done by the service center specified (see section 9.1).  
Never carry out any repairs yourself.

KARL STORZ is liable for safety, reliability and performance of the HF device under the following conditions:

- Full compliance with all instructions regarding the installation and proper use for the intended purpose contained in this instruction manual was maintained.
- Changes, repairs, new settings and similar procedures were carried out only by persons authorized to do this work by KARL STORZ.
- The electrical installations in the relevant room meet the local requirements and statutory provisions.

Fast and satisfactory repairs can only be guaranteed when all required data have been supplied in full.

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The following information is required for returning the device:

- Full address
- Model number
- Serial number
- Software version

▷ Describe the problem, the corresponding application and the accessories used.  
– or –

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### 9 Storage

- ▷ If you store the HF device for longer than one year, pay specific attention to the indicators during automatic functional testing, see section Test for proper functioning, page. 29.
- ▷ Clean the HF device thoroughly before you put it into storage.
- ▷ Store the HF device in a clean, dry place in accordance with the storage conditions.

#### Storage conditions:

- Temperature: -20 °C to +50 °C
- Relative humidity: 0 to 75 %, non-condensing
- Air pressure: 500 to 1060 hPa

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### 9.1 Technical Support

In Germany contact the following service center for maintenance and repair work:

KARL STORZ SE & Co. KG  
Abt. Reparaturservice  
Take-off Gewerbepark 83  
78579 Neuhausen  
Service hotline: +49 7461 708 980  
E-mail: [technicalsupport@karlstorz.com](mailto:technicalsupport@karlstorz.com)

In other countries please contact your local KARL STORZ branch or authorized dealer.

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### 10 Technical data

#### 10.1 Technical data for AUTOCON® III 300

| Insulation type / Classification                                                                       |                                           |
|--------------------------------------------------------------------------------------------------------|-------------------------------------------|
| EMC                                                                                                    | IEC 60601-1-2                             |
| Protection class acc. to IEC 60601-1                                                                   | I                                         |
| Applied part type according to IEC 60601-1                                                             | CF                                        |
| Standard compliance                                                                                    | IEC 60601-1, IEC 60601-1-2, IEC 60601-2-2 |
| Classification according to Annex IX of the EC Directive 93/42/EEC                                     | IIb                                       |
| This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. | <b>CE 0123</b>                            |

**i** The code number after the CE mark indicates the responsible notified body.

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| Programs                         | 300 |
|----------------------------------|-----|
| Number of programs in the device | 300 |
| Individually programmable        | Yes |
| Information shown on the display | Yes |

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| Neutral electrode monitoring                                                                |                                                    |
|---------------------------------------------------------------------------------------------|----------------------------------------------------|
| EASY: Electrode Application System                                                          | Yes                                                |
| Display of one-piece, split and small 'Baby' electrode                                      | Main menu and neutral electrode selection          |
| Contact resistance between individual sections of split neutral electrodes shown on display | Using color and contact indicator                  |
| Lead resistance shown on the display when a non-split neutral electrode is used             | Yes                                                |
| Maximum permissible resistance between the sections of a split electrode                    | 300 Ω                                              |
| Warning signal for hazardous conditions in connection with neutral electrodes               | Visual, acoustic                                   |
| Tones                                                                                       | Warning, activation and button tones, start melody |
| Warning message on the display                                                              | Text message with further information              |

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| Safety features                                               |                                        |
|---------------------------------------------------------------|----------------------------------------|
| ISSys (Integrated Safety System)                              | Yes                                    |
| ARC Control                                                   | ARC CONTROL                            |
| Continuous monitoring of HF leakage current and error message | Text message with further information  |
| Dosage monitoring with error message on the display           | Yes                                    |
| Continuous self-test                                          | Yes                                    |
| Continuous status indication on the display                   | Yes                                    |
| Operating errors shown on the display                         | Text message with further information  |
| System errors shown on the display                            | Text message with further information  |
| Safety inspection                                             | Automatic reminder function (optional) |
| Instruction manual                                            | Hard copy                              |

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| Documentation                                |                                       |
|----------------------------------------------|---------------------------------------|
| Data acquisition and storage in the device   | System information                    |
| Documentation of error states                | Yes                                   |
| Documentation of operating errors            | Yes                                   |
| Retrieval of system messages via the display | Text message with further information |

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| Communication                                                                                 |                            |
|-----------------------------------------------------------------------------------------------|----------------------------|
| Display                                                                                       | Capacitive touch screen 9" |
| USB interface for software updates                                                            | Yes                        |
| External PC interface for service support using KARL STORZ software                           | UART                       |
| External interface for communication between the HF generator and the KARL STORZ OR1™ system. | Ethernet connection        |

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| Service support                                              |     |
|--------------------------------------------------------------|-----|
| Network connection for service support                       | Yes |
| Service support by service programs integrated in the device | Yes |
| Service support via ISSys                                    | Yes |

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| Cooling                    |     |
|----------------------------|-----|
| Convection                 | Yes |
| Temperature-controlled fan | Yes |

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| Operating mode |                                                                                                                                       |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Operating mode | Intermittent 10 s/30 s (on/off)<br><b>i</b> The system switches off automatically following continuous activation lasting 45 seconds. |

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| Specifications                    |                                  |
|-----------------------------------|----------------------------------|
| Max. monopolar power              | 400 W (at 100 Ω)                 |
| Max. bipolar power                | 200 W (at 75 Ω)                  |
| Output frequency                  | 350 kHz / 1 MHz                  |
| Unipolar sockets                  | 1x (footswitch and fingerswitch) |
| Bipolar sockets                   | 2x footswitch                    |
| Connection for footswitch         | 2x                               |
| AUTOSTART                         | Yes                              |
| Scope of supply (see section 3.4) | Inspection protocol              |

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



| Conditions of operation, transport and storage |  | Operation                  | Transport and storage     |
|------------------------------------------------|--|----------------------------|---------------------------|
| Temperature                                    |  | +10 °C to +40 °C           | -20 °C to +60 °C          |
| Relative humidity                              |  | 30 to 75 %, non-condensing | 0 to 90 %, non-condensing |
| Air pressure                                   |  | 700 to 1 060 hPa           | 500 to 1 060 hPa          |
| Operating altitude (max.)                      |  | 3000 m above sea level     |                           |

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




| Mode symbol                                                                           | Designation            | CCS | ARC Control | HF voltage form      | HF power limitation |             | Peak voltage                             | Default values |           |
|---------------------------------------------------------------------------------------|------------------------|-----|-------------|----------------------|---------------------|-------------|------------------------------------------|----------------|-----------|
|                                                                                       |                        |     |             |                      | Effect              | Power range |                                          | Effect         | Max. Watt |
| Unipolar coagulation modes                                                            |                        |     |             |                      |                     |             |                                          |                |           |
|    | Forced Coag<br>1.8 kVp |     |             | sinusoidal modulated | -                   | 1 W – 120 W | 1.8 kVp                                  | ---            | 60        |
|    | Contact Soft           |     |             | sinusoidal constant  | 1<br>2<br>3         | 1 W – 120 W | 250 Vp                                   | 2              | 60        |
|    | Forced coag            |     |             | pulsed modulated     | -                   | 1 W – 80 W  | 3.5 kVp                                  | ---            | 50        |
|   | Resection              |     |             | sinusoidal modulated | -                   | 1 W – 120 W | 2.2 kVp                                  | ---            | 60        |
|  | Spray                  |     |             | pulsed modulated     | 1<br>2<br>3<br>4    | 1 W – 120 W | 3.0 kVp<br>3.8 kVp<br>4.6 kVp<br>5.0 kVp | 2              | 80        |
|  | Forced mixed           |     |             | sinusoidal modulated | 1<br>2<br>3         | 1 W – 120 W | 1.5 kVp<br>2.0 kVp<br>2.5 kVp            | 2              | 60        |
|  | Forced Coag with cut   |     |             | sinusoidal modulated | 1<br>2<br>3<br>4    | 1 W – 250 W | 1.5 kVp<br>1.5 kVp<br>1.3 kVp<br>1.3 kVp | 2              | 80        |

| Mode symbol                                                                         | Designation      | CCS | ARC Control | HF voltage form     | HF power limitation   |                                           | Peak voltage                                   | Default values |           |
|-------------------------------------------------------------------------------------|------------------|-----|-------------|---------------------|-----------------------|-------------------------------------------|------------------------------------------------|----------------|-----------|
|                                                                                     |                  |     |             |                     | Effect                | Power range                               |                                                | Effect         | Max. Watt |
| Bipolar cutting modes                                                               |                  |     |             |                     |                       |                                           |                                                |                |           |
|  | Bip. cutting     | Yes | Yes         | sinusoidal constant | -                     | 1 W – 200 W                               | 400 Vp                                         | ---            | 100       |
|  | Bipolar scissors |     |             | sinusoidal constant | -                     | 1 W – 120 W                               | 200 Vp                                         | ---            | 40        |
|  | VAP CUT          | Yes | Yes         | sinusoidal constant | 1<br>2<br>3<br>4<br>5 | 120 W<br>200 W<br>200 W<br>200 W<br>200 W | 290 Vp<br>330 Vp<br>380 Vp<br>440 Vp<br>500 Vp | 3              | ---       |

| Mode symbol                                                                         | Designation   | CCS | ARC Control | HF voltage form     | HF power limitation |                      | Peak voltage               | Default values |           |
|-------------------------------------------------------------------------------------|---------------|-----|-------------|---------------------|---------------------|----------------------|----------------------------|----------------|-----------|
|                                                                                     |               |     |             |                     | Effect              | Power range          |                            | Effect         | Max. Watt |
| Bipolar coagulation modes                                                           |               |     |             |                     |                     |                      |                            |                |           |
|  | RoBi®         |     |             | sinusoidal constant | -                   | 1 W – 100 W          | 110 Vp                     | ---            | 40        |
|  | Standard plus |     |             | sinusoidal constant | -                   | 1 W – 120 W          | 150 Vp                     | ---            | 50        |
|  | Standard      |     |             | sinusoidal constant | -                   | 1 W – 120 W          | 150 Vp                     | ---            | 40        |
|  | VAP COAG      |     |             | sinusoidal constant | 1<br>2<br>3         | 40 W<br>60 W<br>80 W | 110 Vp<br>150 Vp<br>550 Vp | 2              | ---       |




| Mode symbol                                                                          | Designation      | CCS | ARC Control | HF voltage form | HF power limitation  |             | Peak voltage | Default values |           |
|--------------------------------------------------------------------------------------|------------------|-----|-------------|-----------------|----------------------|-------------|--------------|----------------|-----------|
|                                                                                      |                  |     |             |                 | Effect               | Power range |              | Effect         | Max. Watt |
| Bipolar coagulation modes                                                            |                  |     |             |                 |                      |             |              |                |           |
|   | Standard AUTO    |     |             |                 | sinusoidal constant  | -           | 5 W – 120 W  | 150 Vp         | ---<br>40 |
|   | Bipolar scissors |     |             |                 | sinusoidal constant  | -           | 1 W – 120 W  | 200 Vp         | ---<br>40 |
|   | Micro            |     |             |                 | sinusoidal constant  | -           | 0.1 W – 40 W | 90 Vp          | ---<br>10 |
|  | Forced           |     |             |                 | sinusoidal modulated | -           | 1 W – 100 W  | 550 Vp         | ---<br>50 |

**i** These max. values are not necessarily created at rated load.  
The HF power limitation is subject to a tolerance of  $\pm 20\%$ .

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

### 10.2 Output, voltage and current diagrams Unipolar cutting – Cut reduced 200W

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Cut reduced 200 W' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Cut reduced 200 W'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Cut reduced 200 W' = 200 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Cut reduced 200 W' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 400    |
| 2      | 450    |
| 3      | 560    |
| 4      | 650    |
| 5      | 650    |
| 6      | 700    |
| 7      | 700    |
| 8      | 700    |
| 9      | 750    |

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

#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Cut' = 150 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Cut'. Rated load resistance = 200 Ω

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**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Cut' = 300 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Cut' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 400    |
| 2      | 450    |
| 3      | 560    |
| 4      | 650    |
| 5      | 650    |
| 6      | 700    |
| 7      | 700    |
| 8      | 700    |
| 9      | 750    |

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

#### Measurement at ohmic resistances

Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Cut reduced 50 W' = 25 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Cut reduced 50 W'. Rated load resistance = 500 Ω

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**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Cut reduced 50 W' = 50 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Cut reduced 50W' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 280    |
| 2      | 340    |
| 3      | 380    |
| 4      | 400    |
| 5      | 400    |
| 6      | 400    |
| 7      | 450    |
| 8      | 450    |
| 9      | 450    |

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Resection'

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- Table of power output P [W] as a function of the setting 'Unipolar Cutting Resection'. Rated load resistance = 500 Ω

| Efecto | P (W) |
|--------|-------|
| 1      | 250   |
| 2      | 250   |
| 3      | 250   |
| 4      | 250   |
| 5      | 250   |

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
- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Resection' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 650    |
| 2      | 700    |
| 3      | 700    |
| 4      | 700    |
| 5      | 750    |

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Cut 300-400 W'

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- Table of power output P [W] as a function of the setting 'Unipolar Cutting Cut 300-400 W'. Rated load resistance = 100 Ω

| Efecto | P (W) |
|--------|-------|
| 1      | 300   |
| 2      | 350   |
| 3      | 400   |

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Cut 300-400 W' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 650    |
| 2      | 650    |
| 3      | 650    |

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

#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Cut blend' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Cutting Cut blend'. Rated load resistance = 500 Ω

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**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Cut blend' = 200 W

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
- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Cut blend' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 1400   |
| 2      | 1400   |
| 3      | 1400   |
| 4      | 1400   |
| 5      | 1500   |
| 6      | 1600   |
| 7      | 1600   |
| 8      | 1600   |
| 9      | 1600   |

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 400 W 1',  
coag phase

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
- Table of power output P [W] as a function of the setting 'Unipolar Cutting Intermittent 400 W 1'.  
Rated load resistance = 500 Ω

| Efecto | P (W) Coag | P (W) Fase Cut |
|--------|------------|----------------|
| 1      | 17         | 400            |
| 2      | 21         | 400            |
| 3      | 26         | 400            |
| 4      | 30         | 400            |
| 5      | 35         | 400            |

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 400W 1',  
cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Intermittent 400W 1' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 750    |
| 2      | 750    |
| 3      | 750    |
| 4      | 750    |
| 5      | 750    |

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 400 W 2'  
coag phase

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
- Table of power output P [W] as a function of the setting 'Unipolar Cutting Intermittent 400 W 2'.  
Rated load resistance = 500 Ω

| Efecto | P (W) Coag | P (W) Fase Cut |
|--------|------------|----------------|
| 1      | 17         | 400            |
| 2      | 21         | 400            |
| 3      | 26         | 400            |
| 4      | 30         | 400            |
| 5      | 35         | 400            |

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**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 400 W 2'  
cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Intermittent 400 W 2' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 750    |
| 2      | 750    |
| 3      | 750    |
| 4      | 750    |
| 5      | 750    |

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 400W 3'  
coag phase

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
- Table of power output P [W] as a function of the setting 'Unipolar Cutting Intermittent 400W 3'. Rated load resistance = 500 Ω

| Efecto | P (W) Coag | P (W) Fase Cut |
|--------|------------|----------------|
| 1      | 17         | 400            |
| 2      | 21         | 400            |
| 3      | 26         | 400            |
| 4      | 30         | 400            |
| 5      | 35         | 400            |

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**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 400 W 3'  
cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Intermittent 400 W 3' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 750    |
| 2      | 750    |
| 3      | 750    |
| 4      | 750    |
| 5      | 750    |

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [ $\Omega$ ] for the setting 'Unipolar Cutting Intermittent 300W 1',  
coag phase

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
- Table of output power P [W] as a function of the setting 'Unipolar Cutting Intermittent 300W 1'.  
Rated load resistance = 500  $\Omega$

| Efecto | P (W) Coag | P (W) Fase Cut |
|--------|------------|----------------|
| 1      | 17         | 300            |
| 2      | 21         | 300            |
| 3      | 26         | 300            |
| 4      | 32         | 300            |
| 5      | 40         | 300            |

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**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 300 W 1'  
cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Intermittent 300 W 1' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 650    |
| 2      | 650    |
| 3      | 650    |
| 4      | 700    |
| 5      | 750    |

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
**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 300 W 2' coag phase

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
- Table of output power P [W] as a function of the setting 'Unipolar Cutting Intermittent 300 W 2'.  
Rated load resistance = 500 Ω

| Efecto | P (W) Coag | P (W) Fase Cut |
|--------|------------|----------------|
| 1      | 17         | 300            |
| 2      | 21         | 300            |
| 3      | 26         | 300            |
| 4      | 32         | 300            |
| 5      | 40         | 300            |

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**Measurement at ohmic resistances**

Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 300 W 2' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Intermittent 300 W 2' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 650    |
| 2      | 650    |
| 3      | 650    |
| 4      | 700    |
| 5      | 750    |

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [ $\Omega$ ] for the setting 'Unipolar Cutting Intermittent 300 W 3'  
coag phase

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
- Table of output power P [W] as a function of the setting 'Unipolar Cutting Intermittent 300 W 3'.  
Rated load resistance = 500  $\Omega$

| Efecto | P (W) Coag | P (W) Fase Cut |
|--------|------------|----------------|
| 1      | 17         | 300            |
| 2      | 21         | 300            |
| 3      | 26         | 300            |
| 4      | 32         | 300            |
| 5      | 40         | 300            |

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Cutting Intermittent 300 W 3' cut phase

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Cutting Intermittent 300 W 3' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 650    |
| 2      | 650    |
| 3      | 650    |
| 4      | 700    |
| 5      | 750    |

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
Unipolar coagulation – Forced Coag 1.8 kVp

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced Coag 1.8 kVp'  
= 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced Coag 1.8 kVp'.  
Rated load resistance = 500 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced Coag 1.8 kVp'  
= 120 W

- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Forced Coag 1.8 kVp' (idle mode) = 1800 Vp

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
Unipolare Koagulation – Kontakt Soft

Unipolar coagulation – Contact Soft

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [ $\Omega$ ] for the setting 'Unipolar Coagulation Contact Soft' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Contact Soft'.  
Rated load resistance = 75  $\Omega$

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Contact Soft' = 120 W

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
- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Contact Soft' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 250    |
| 2      | 250    |
| 3      | 250    |

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced Coag' = 40 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced Coag'.  
Rated load resistance = 1000 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced Coag' = 80 W

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- HF output voltage U [Vp] with the setting 'Unipolar Coagulation Forced Coag' (idle mode) = 3500 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Resection' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Resection'. Rated load resistance = 500 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Resection' = 120 W

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

#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Spray' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Spray'. Rated load resistance = 500 Ω

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Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Spray' = 120 W

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- Diagram of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Spray' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 3000   |
| 2      | 3800   |
| 3      | 4600   |
| 4      | 5000   |

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
Unipolare Koagulation – Forciert mixed

Unipolar coagulation – Forced mixed

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [ $\Omega$ ] for the setting 'Unipolar Coagulation Forced mixed' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced mixed'. Rated load resistance = 500  $\Omega$

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced mixed' = 120 W

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- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Forced mixed' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 1500   |
| 2      | 2000   |
| 3      | 2500   |

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
Unipolar coagulation – Forced Coag with Cut

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
- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Unipolar Coagulation Forced Coag with Cut'  
= 125 W

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
- Diagram of power output P [W] as a function of the setting 'Unipolar Coagulation Forced Coag with Cut'.  
Rated load resistance = 500 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [ $\Omega$ ] for the setting 'Unipolar Coagulation Forced Coag with Cut'  
= 250 W

| Effect | U (Vp) |
|--------|--------|
| 1      | 1500   |
| 2      | 1500   |
| 3      | 1300   |
| 4      | 1300   |

- Table of HF output voltage U [Vp] as a function of the setting 'Unipolar Coagulation Forced Coag with Cut' (idle mode)

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

#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Cutting' = 100 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Cutting Bip. Cutting'.  
Rated load resistance = 75 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. Cutting' = 200 W

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- HF output voltage U [Vp] with the setting 'Bipolar Cutting Bip. Cutting' (idle mode) = 400 Vp

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
Corte bipolar – Tijeras bipolar.

Measurement at ohmic resistances


- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bipolar Scissors' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Cutting Bipolar Scissors'. Rated load resistance = 75 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bipolar Scissors' = 120 W

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- HF output voltage U [Vp] with the setting 'Bipolar Cutting Bipolar Scissors' (idle mode) = 200 Vp

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
**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Cutting Bip. VAP CUT'

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- Table of power output P [W] as a function of the setting 'Bipolar Cutting VAP CUT'. Rated load resistance = 75 Ω

| Efecto | P (W) |
|--------|-------|
| 1      | 120   |
| 2      | 200   |
| 3      | 200   |
| 4      | 200   |
| 5      | 200   |

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- Table of HF output voltage U [Vp] as a function of the setting "Bipolar Cutting VAP CUT" (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 290    |
| 2      | 330    |
| 3      | 380    |
| 4      | 440    |
| 5      | 500    |

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

**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation RoBi®' = 50 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation RoBi®'. Rated load resistance = 25 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation RoBi®' = 100 W

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- HF output voltage U [Vp] with the setting 'Bipolar Coagulation RoBi®' (idle mode) = 110 Vp

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
Bipolare Koagulation – Standard plus

Bipolar coagulation – Standard plus

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
Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard plus' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Standard plus'. Rated load resistance = 50 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard plus' = 120 W

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- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Standard plus' (idle mode) = 150 Vp

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

**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Standard'. Rated load resistance = 50 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard' = 120 W

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- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Standard' (idle mode) = 150 Vp

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation VAP COAG'

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- Table of power output P [W] as a function of the setting 'Bipolar Coagulation VAP COAG'. Rated load resistance = 50 Ω (at effect 1 and 2)  
Rated load resistance = 100 Ω (at effect 3)

| Efecto | P (W) |
|--------|-------|
| 1      | 40    |
| 2      | 60    |
| 3      | 80    |

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- Table of HF output voltage U [Vp] as a function of the setting 'Bipolar Coagulation VAP COAG' (idle mode)

| Efecto | U (Vp) |
|--------|--------|
| 1      | 110    |
| 2      | 150    |
| 3      | 550    |

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
Bipolare Koagulation – Standard AUTO

Bipolar coagulation – Standard AUTO

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard AUTO' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Standard AUTO'. Rated load resistance = 50 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Standard AUTO' = 120 W

- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Standard AUTO' (idle mode) = 150 Vp

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
Bipolare Koagulation – Bipolare Schere

Bipolar coagulation – Bipolar scissors

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
#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bipolar Scissors' = 60 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Bipolar Scissors'. Rated load resistance = 75 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Bipolar Scissors' = 120 W

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- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Bipolar Scissors' (idle mode) = 200 Vp

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

Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Micro' = 20 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Micro'. Rated load resistance = 50 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Micro' = 40 W

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- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Micro' (idle mode) = 90 Vp

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

**Measurement at ohmic resistances**

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Forced' = 50 W

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
- Diagram of power output P [W] as a function of the setting 'Bipolar Coagulation Forced'. Rated load resistance = 50 Ω

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#### Measurement at ohmic resistances

- Diagram of power output P [W] as a function of the load resistance R [Ω] for the setting 'Bipolar Coagulation Forced' = 100 W

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- HF output voltage U [Vp] with the setting 'Bipolar Coagulation Forced' (idle mode) = 550 Vp

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### 10.3 Technical documentation

On request, the manufacturer will provide those service manuals, itemized parts listings, descriptions, sets of adjustment instructions and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications, in particular in the context of ongoing development and improvement of our devices.

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## 11 Accessories/Replacement parts

Original KARL STORZ accessories are suitable for use with the AUTOCON® III 300 devices. When using accessories made by other manufacturers, the user must ensure that they are designed for and compatible with the maximum HF peak voltage of the HF device.

For the use and correct reprocessing of the autoclavable devices, compliance with the relevant instruction manuals accompanying these devices is required.

Detailed information on accessories and replacement parts is available in the corresponding KARL STORZ catalogs.

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## 12 EMC

The AUTOCON® III 300 is suitable for use in Professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

### 12.1 Guidance and manufacturer's declaration as per IEC 60601-1-2

**WARNING:** Medical electrical devices are subject to special precautions regarding electromagnetic compatibility (EMC). Observe the EMC instructions in this section during installation and commissioning.

**WARNING:** The AUTOCON® III 300 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device or system should be observed to verify normal operation in the configuration in which it will be used.

**WARNING:** The use of portable and mobile HF communication equipment may have an impact on this or other medical electrical devices. Such devices including their accessories such as e.g. antennas should not be operated closer than 30 cm (12 inches) to the AUTOCON® III 300 and attached accessories. Otherwise, degradation of the performance of this equipment could result.

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**WARNING:** The use of an accessory or cable with the AUTOCON® III 300 other than those specified in this manual may result in increased emissions or decreased immunity of the AUTOCON® III 300. When using an accessory or cable with AUTOCON® III 300 other than those specified in this manual (section 11), it becomes the responsibility of the user of the AUTOCON® III 300 to determine compliance with IEC 60601-1-2 when using this item.

The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

The tables and guidelines that are included in this section provide information to the customer or user that is essential in determining the suitability of the device or system for the electromagnetic environment of use, and in managing the electromagnetic system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If electromagnetic interference arises while using the device, the user can eliminate it by taking the following actions:

- ▷ reorient or relocate the receiving device
- ▷ increase the separation between the equipment
- ▷ connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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### Essential performance

The standard IEC 60601-2-2 covers the requirements for basic safety of HF surgical devices in accordance with IEC 60601-1.

As part of the risk management process, the following additional essential performance characteristics were identified:

- The accuracy of the maximum HF output voltage and the effect on patient tissue, which are summarized in the specific standard 60601-2-2.
- A stable neutral electrode monitoring system.
- The minimum noise level of the activation indicator and the acoustic notification sounds.

Due to the EMC ambient conditions beyond the levels listed in the tables below, the function of the essential performance may be impaired or fail.

This may cause malfunction or failure of the neutral electrode monitoring system as well as the acoustic activation indications and notification sounds. The accuracy of the RF output power may deviate more than the permissible 20 %.

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| Guidance and manufacturer's declaration – electromagnetic emissions (IEC 60601-1-2)                                                                                                          |            |                                                                                                                                                                                                                                                                        |  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| AUTOCON® III 300 is intended for use in the electromagnetic environment specified below. The customer or user of the AUTOCON® III 300 must make sure that it is used in such an environment. |            |                                                                                                                                                                                                                                                                        |  |
| Interference emissions measurements                                                                                                                                                          | Compliance | Electromagnetic environment – guidelines                                                                                                                                                                                                                               |  |
| HF emissions acc. to CISPR 11                                                                                                                                                                | Group 2    | The AUTOCON® III 300 must emit electromagnetic energy in order to perform its intended function. Nearby electronic devices may be affected.                                                                                                                            |  |
| HF emissions acc. to CISPR 11                                                                                                                                                                | Class A    | The AUTOCON® III 300 is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.<br>– Professional healthcare facility environment – |  |
| Harmonic emissions acc. to IEC 61000-3-2                                                                                                                                                     | Class A    |                                                                                                                                                                                                                                                                        |  |
| Voltage fluctuations/flicker emissions acc. to IEC 61000-3-3                                                                                                                                 | Complies   |                                                                                                                                                                                                                                                                        |  |

| Table 4 – ENCLOSURE PORT (according to IEC 60601-1-2)      |                                   |                                                  |                                                                                     |
|------------------------------------------------------------|-----------------------------------|--------------------------------------------------|-------------------------------------------------------------------------------------|
| Phenomenon                                                 | Basic EMC standard or test method | IMMUNITY TEST LEVELS                             |                                                                                     |
|                                                            |                                   | Professional healthcare facility environment     | HOME HEALTHCARE ENVIRONMENT                                                         |
| ELECTROSTATIC DISCHARGE                                    | IEC 61000-4-2                     | ±8 kV contact<br>±2 kV, ±4 kV, ±8 kV, ±15 kV air |                                                                                     |
| Radiated RF EM fields                                      | IEC 61000-4-3                     | 3 V/m<br>80 MHz – 2.7 GHz<br>80 % AM at 1 kHz    | Not applicable since released only for professional healthcare facility environment |
| Proximity fields from RF wireless communications equipment | IEC 61000-4-3                     | See table 9                                      |                                                                                     |
| RATED power frequency magnetic fields                      | IEC 61000-4-8                     | 30 A/m<br>50 Hz or 60 Hz                         |                                                                                     |

**Table 5 – Input a.c. power PORT** (according to IEC 60601-1-2)

| Phenomenon                                  | Basic EMC standard | IMMUNITY TEST LEVELS                                                                                                                                              |                                                                                     |
|---------------------------------------------|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                                             |                    | Professional healthcare facility environment                                                                                                                      | HOME HEALTHCARE ENVIRONMENT                                                         |
| Electrical fast transients/bursts           | IEC 61000-4-4      | ±2 kV<br>100 kHz repetition frequency                                                                                                                             |                                                                                     |
| Surges<br>Line-to-line                      | IEC 61000-4-5      | ±0.5 kV, ±1 kV                                                                                                                                                    |                                                                                     |
| Surges<br>Line-to-ground                    | IEC 61000-4-5      | ±0.5 kV, ±1 kV, ±2 kV                                                                                                                                             |                                                                                     |
| Conducted disturbances induced by RF fields | IEC 61000-4-6      | 3 V<br>0.15 MHz – 80 MHz<br>6 V in ISM bands between 0.15 MHz and 80 MHz<br>80 % AM at 1 kHz                                                                      | Not applicable since released only for professional healthcare facility environment |
| Voltage dips                                | IEC 61000-4-11     | 0 % $U_T$ ; 0.5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br><br>0 % $U_T$ ; 1 cycle and<br>70 % $U_T$ ; 25/30 cycles<br>Single phase: at 0° Grad |                                                                                     |
| Voltage interruptions                       | IEC 61000-4-11     | 0 % $U_T$ ; 250/300 cycles                                                                                                                                        |                                                                                     |

**Table 8 – Signal input/output parts PORT** (according to IEC 60601 -1-2)

| Phenomenon                                  | Basic EMC standard | IMMUNITY TEST LEVELS                                                                         |                                                                                     |
|---------------------------------------------|--------------------|----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                                             |                    | Professional healthcare facility environment                                                 | HOME HEALTHCARE ENVIRONMENT                                                         |
| ELECTROSTATIC DISCHARGE                     | IEC 61000-4-2      | ±8 kV contact<br>±2 kV, ±4 kV, ±8 kV, ±15 kV air                                             |                                                                                     |
| Electrical fast transients/bursts           | IEC 61000-4-4      | ±1 kV<br>100 kHz repetition frequency                                                        |                                                                                     |
| Surges<br>Line-to-ground                    | IEC 61000-4-5      | ±2 kV                                                                                        |                                                                                     |
| Conducted disturbances induced by RF fields | IEC 61000-4-6      | 3 V<br>0.15 MHz – 80 MHz<br>6 V in ISM bands between 0.15 MHz and 80 MHz<br>80 % AM at 1 kHz | Not applicable since released only for professional healthcare facility environment |

| Table 9 – Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment (according to IEC 60601-1-2)                                                                           |               |                                                                                |                                       |                      |                 |                              |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|--------------------------------------------------------------------------------|---------------------------------------|----------------------|-----------------|------------------------------|
| Test frequency<br>(MHz)                                                                                                                                                                                  | Band<br>(MHz) | Service                                                                        | Modulation                            | Maximum Power<br>(W) | Distance<br>(m) | IMMUNITY TEST LEVEL<br>(V/m) |
| 385                                                                                                                                                                                                      | 380 – 390     | TETRA 400                                                                      | Pulse modulation<br>18 Hz             | 1.8                  | 0.3             | 27                           |
| 450                                                                                                                                                                                                      | 430 – 470     | GMRS 460,<br>FRS 460                                                           | FM<br>± 5 kHz deviation<br>1 kHz sine | 2                    | 0.3             | 28                           |
| 710                                                                                                                                                                                                      | 704 – 787     | LTE Band 13, 17                                                                | Pulse modulation<br>217 Hz            | 0.2                  | 0.3             | 9                            |
| 745                                                                                                                                                                                                      |               |                                                                                |                                       |                      |                 |                              |
| 780                                                                                                                                                                                                      |               |                                                                                |                                       |                      |                 |                              |
| 810                                                                                                                                                                                                      | 800 – 960     | GSM 800/900,<br>TETRA 800,<br>iDEN 820,<br>CDMA 850,<br>LTE Band 5             | Pulse modulation<br>18 Hz             | 2                    | 0.3             | 28                           |
| 870                                                                                                                                                                                                      |               |                                                                                |                                       |                      |                 |                              |
| 930                                                                                                                                                                                                      |               |                                                                                |                                       |                      |                 |                              |
| 1 720                                                                                                                                                                                                    | 1 700 – 1 990 | GSM 1800;<br>CDMA 1900;<br>GSM 1900;<br>DECT;<br>LTE Band 1, 3,<br>4, 25; UMTS | Pulse modulation<br>217 Hz            | 2                    | 0.3             | 28                           |
| 1 845                                                                                                                                                                                                    |               |                                                                                |                                       |                      |                 |                              |
| 1 970                                                                                                                                                                                                    |               |                                                                                |                                       |                      |                 |                              |
| 2 450                                                                                                                                                                                                    | 2 400 – 2 570 | Bluetooth,<br>WLAN 802.11 b/g/n,<br>RFID 2450,<br>LTE Band 7                   | Pulse modulation<br>217 Hz            | 2                    | 0.3             | 28                           |
| 5 240                                                                                                                                                                                                    | 5 100 – 5 800 | WLAN 802.11 a/n                                                                | Pulse modulation<br>217 Hz            | 0.2                  | 0.3             | 9                            |
| 5 500                                                                                                                                                                                                    |               |                                                                                |                                       |                      |                 |                              |
| 5 785                                                                                                                                                                                                    |               |                                                                                |                                       |                      |                 |                              |
| NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the AUTOCON® III 300 may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3. |               |                                                                                |                                       |                      |                 |                              |

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### 13 Disposal

**i** The national regulations must be observed when disposing of or recycling the device or its components.

| Symbol                                                                              | Designation                                                                                                                                                           |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Devices bearing this symbol must be collected separately as electrical and electronic devices. Within the EU, disposal is handled free of charge by the manufacturer. |

▷ For questions on the disposal of the device, please contact the Technical Support; see section Technical Support on page 68.

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7171 Millcreek Drive, Mississauga, Ontario

LSN 3R3 Canada

Phone: +1 905 816-4500, Fax: +1 905 816-4599  
Toll free phone: 1-800-268-4880 (Canada only)  
Toll free fax: 1-800-482-4198 (Canada only)  
E-Mail: info-canada@karlstorz.com

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E-Mail: info-az@karlstorz.com

KARL STORZ ENDOSKOPE –

East Mediterranean and Gulf (Offshore) S.A.L.

Spark Tower 1<sup>st</sup> floor

Charlous Helou St, Horch Tabet – Sin Efil

Beirut, Lebanon

Phone: +961 1 501105, Fax: +961 1 501950

E-Mail: info@karlstorz-emg.com

**Sociedades distribuidoras**

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch  
Room 803-805, 8F, Jin Jiang International Building  
No. 1 West Linjiang Road, Wuhou District,  
6100414, Chengdu, People's Republic of China  
Phone: +86 28 86587977, Fax: +86 28 86587975  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch  
Room 2001-2005, 20F N-MEDIA International Center,  
No. 167 Youth Avenue, Shenhe District,  
110014, Shenyang, People's Republic of China  
Phone: +86 24 23181118, Fax: +86 24 23181119  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch  
Room 02B & 03 & 04A, 35F Teem Tower,  
No. 208 Tianhe Road, Tianhe District,  
510620, Guangzhou, People's Republic of China  
Phone: +86 20 87321281, Fax: +86 20 87321286  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd  
No. 8 Commonwealth Lane #03-02  
Singapore 149555, Singapore  
Phone: +65 69229150, Fax: +65 69229155  
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.  
9F Hyowon-Building  
97, Jungdae-ro, Songpa-gu  
05719 Seoul, Korea  
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299  
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.  
12F, No. 192, Sec. 2, Chung Hsin Rd.,  
Sindian District, New Taipei City, Taiwan  
Phone: +886 933 014 160, Fax: +886 2 8672 6399  
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG  
Representative Office Philippines  
1901 Picadilly Star Bldg., 4<sup>th</sup> Avenue, BGC  
Taguig City 1636, Philippines  
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11  
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.  
Stage Bldg. 8F, 2-7-2 Fujimi  
Chiyoda-ku, Tokyo 102-0071, Japan  
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633  
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.  
68 Waterloo Road, Macquarie Park NSW 2113  
P O Box 50 Lane Cove NSW 1595, Australia  
Phone: +61 (0)2 9490 6700  
Toll free: 1800 996 562 (Australia only)  
Fax: +61 (0)2 9420 0695  
E-Mail: karlstorz@karlstorz.com.au  
www.karlstorz.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.  
P.O. 6061, Roggebaai 18012  
Cape Town, South Africa  
Phone: +27 21 417 2600, Fax: +27 21 421 5103  
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasakhstan  
Saykara, 6, BC "Arman", off. 910  
010000 Astana, Republic of Kazakhstan  
Phone: +7 7172 552-549, 552-788, Fax: -444  
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE  
East Mediterranean & Gulf (branch)  
Building West Side 7A – Unit 7MA – 3008  
Dubai Airport Free Zone, P.O. Box 54983  
Dubai - United Arab Emirates  
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282  
Service Hotline: +971 (0)4 3415882  
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited  
11<sup>th</sup> Floor, Dr. Gopal Das Bhawan  
28, Barakhamba Road  
New Delhi 110001, India  
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010  
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG  
Interchange 21 Tower, Level 33  
399 Sukhumvit Road,  
North Klongtoey, Wattana,  
10110 Bangkok, Thailand  
Phone: +84 28 3823 8000  
Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG  
Resident Representative Office  
14<sup>th</sup> Floor, MPHase Saigon  
39 Le Duan, District 1  
Ho Chi Minh City, Vietnam  
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039  
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.  
Room 2503-05, 25F AXA Tower, Landmark East,  
No. 100 How Ming Street, Kwon Tong, Kowloon,  
Hong Kong, People's Republic of China  
Phone: +852 28 65 2411, Fax: +852 28 65 4114  
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch  
Room 1805-1807, Building B, 18F Beijing IFC  
No. 8, Jiaquomenwai Street, Chaoyang District,  
100022, Beijing, People's Republic of China  
Phone: +86 10 5638188, Fax: +86 10 5638199  
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch  
Room 701A Building 5 & Room 501 Building 7,  
No. 3000 Longgong Avenue, Pilot Free Trade Zone,  
201203, Shanghai, People's Republic of China  
Phone: +86 21 60339888, Fax: +86 21 60339808  
E-Mail: info@karlstorz.com.cn

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