FOR INFORMATION ONLY.
WHEN PERFORMING
THE ASSAY ALWAYS REFER
TO PACKAGE INSERT
SUPPLIED
WITH THE KIT

RO

CanAg CA125 EIA

REF

400-10

IVD

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EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIIFICAȚIA SIMBOLURILOR
RU	ОБОНАЧЕНИЯ
SV	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/
Holdbar til/Verwendbar bis/
Ημερομηνία λήξης/Fecha
de caducidad/Kölblük kuni/
Utiliser jusque/Rok valjanosti/
Felhasználható/Utilizzare entro/
Sunautoti iki/Izletot līdz/Houdbaar
tot/Brukes innen/Użyty przed/
Prazo de validade/Expiră la/
Использовать до/Använd före/
Použite ne do/ Uporabno do/
Upotrebljivo do/Son Kullanma Tarihi

LOT

Batch code/Номер на партида/
Číslo šarže/Lotnummer/
Chargenbezeichnung/Αριθμός
Παρτίδας/Código de lote/Partii
kood/Code du lot/Kod serije/
Sarzszzám/Codice del lotto/
Partijos kodas/Partijas kods/Lot
nummer/Partikode/Kod partii/
Código do lote/Număr de lot/
Номер лота/Lotnummer/Číslo
šarže/Številka serije/Kod partije/
Parti Kodu

Date of manufacture/Dاتا на производство/Datum výroby/
Produktionsdato/Herstellungsdatum/
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/
Date de fabrication/Datum proizvodnje/
Gyártási idő/Data di produzione/
Pagaminimo data/Ražošanas datums/
Productiedatum/Fremstillingsdato/
Data produkcji/Data de fabrico/Data fabricației/Дата производства/
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Üretim tarihi

Temperature limitation/
Температурни граници/
Teplotní omezení/
Temperaturbegrænsning/
Temperaturbegrenzung/
Περιορισμοί θερμοκρασίας/
Limites de temperatura/
Temperatuuri piirang/
Limite de température/
Temperaturno ograničenje/
Hőmérsékletre vonatkozó korlátozás/
Limiti di temperatura/
Temperatūriniai apribojimai/
Temperatūras ierobežojums/
Temperaturbepæring/
Temperaturbegrensinger/
Temperatury granicne/
Limite de temperatura/
Limite de temperatură/
Температурный режим/
Temperaturbegrænsning/
Teplotné obmedzenie
Omejitve temperature/
Temperaturno ograničenje/
Sıcaklık sınırlaması/

IVD

In Vitro Diagnostic Medical Device/
Медицински уред за диагностика
ин vitro/Diagnostický zdravotnický
prostředek in vitro/Medicinsk udstyr til
in vitro-diagnostik/In-vitro-Diagnostikum/
Ιατροτεχνολογικό προϊόν για διάγνωση
In Vitro/Dispositivo médico para
diagnóstico in vitro/In vitro diagnostiline
meditsiiniseade/Dispositif médical de
diagnostic in vitro/Diagnostički medicinski
uređaj In Vitro/In vitro orvodisagnostikai
eszköz/Dispositivo medico per test
diagnostici in vitro/In Vitro Diagnostiné
Medicinos Priemonė/Mediciniska ierice
in vitro diagnostikai/In vitro-diagnostisch
medisch instrument/In vitro diagnostisk
medisinsk utstyr/Wyrób medyczny do
diagnostyki in vitro/Dispositivo Médico
de Diagnóstico In Vitro/Dispozitiv medical
pentru diagnostic in vitro/Только для
диагностики In Vitro/Endast för in
vitro-diagnostik/ Zdravotnicka pomůcka na
diagnostiku in vitro/In vitro diagnostični
pripomoček/Diagnostički medicinski
uređaj In Vitro/<96> testleri için yeterlilik
içerir

Contains sufficient for <96> tests/Съдържа
достатъчно количество за тестове
<96>/Lze použít pro <96> testů/Ineholder
tilstrækkeligt/Inhalt ausreichend für <96>
Prüfungen/Περιεχόμενο επαρκές για
«96» εξετάσεις/Contenido suficiente para
<96> ensayos/Kogusest piisab <96> testi
läbiviimiseks/Contenu suffisant pour "96"
tests/Sadržj dovoljno za <96> testova/A
doboz tartalma <96> vizsgálat elvégzéséhez
elegendő/Contenuto sufficiente per "96"
saggi/Turinys skirtas atlikti <96> tyrimus/
Saturis pietiekams <96> testiem/Inhoud
voldoende voor "96" testen/til "96" test/
Tilstrækkelig innhold for <96> prøver/
Wystarczy na wykonanie <96> testów/
Conteúdo suficiente para "96" ensaios/
Conținut suficient pentru 96 de teste/
Содержит достаточные количества для
«96» определений/Innehåller tillräckligt
till "96" antal tester/Obsah postačuje na
tento počet testov: <96>/Vsebinsa zadostuje
za <96> testov/Sadržina dovoljna za <96>
testova/<96> testleri için yeterlilik içerir

REF

Catalogue number/Каталожен номер/
Katalogové číslo/Katalognummer/
Bestellnummer/Αριθμός καταλόγου/
Número de catálogo/Katalogoi number/
Numéro de catalogue/Kataloški broj/
Katalógusszám/Numero di catalogo/
Katalogo numeris/Numurs katalogā/
Catalogusnummer/Katalognummer/
Numer katalogowy/Número do catálogo/
Număr de catalog/Номер по каталогу/
Produktnummer/Katalogové číslo/
Kataloška številka/Kataloški broj/
Katalog numarası

Consult Instructions for Use/
Прочетете инструкцията за
употреба/Konzultujte s návodem
k použití/Se brugsanvisning/Siehe
Gebrauchsanweisung/Συμβουλευτείτε
πες Οδηγίες σχετικά με τη χρήση/
Consulte las instrucciones de uso/
Vt kasutusjuhendit/Consulter le mode
d'emploi/Pročítajte upute za uporabu/
Olvassa el a használati utasítást/
Consultare le istruzioni per l'uso/Dél
naudojimo žiūrėkite instrukcijas/Izlasiet
lietošanas instrukciju/Raadpleeg de
instructies voor gebruik/Les instruksene
fer bruk/Sprawdzić w instrukcji użycia/
Consulte as Instruções de Utilização/
Consultati instrucțiunile de utilizare/
Обратитесь к инструкции по
применению/Se brugsanvisning/
Prečítajte si návod na používanie/
Pročítajte uputstvo za upotrebu/
Kullanım Talimatlarına Bakınız

CONT

Contents of kit/Съдържание на набора/
Obsah soupravy/Kittets indhold/Inhalt
des Kits/Περιεχόμενα του kit/Contenido
del kit/Komplekt sisaldab/Contenu du
kit/Sadržaj opreme/A készlet tartalma/
Contenuto del kit/Rinkinio turinys/
Komplekta saturs/Inhoud van de set/
Settets innhold/Zawartość zestawu/
Conteúdo do kit/Conținutul setului/
Компоненты набора/Kit innehåll/
Obsah súpravy/Vseбина kompleta/Sadržaj
opreme/Kitin içindekiler

Biological risks/Биологическа
опасност/Biológická rizika/Biologisk
fare/Biologische Gefahren/Βιολογικοί
κίνδυνοι/Riesgos biológicos/
Biologilised ohud/Risques biologiques/
Biolóskli rizici/Biológiai kockázatok/Rischi
biologici/Biologinis pavojus/Biologiskais
risks/Biologische risico's/Biologische
risikoer/Zagrozenie biologiczne/Riscos
biológicos/ Biologisk risk/Pericole
biologice/Биологическая опасность/
Biologický rizikové/Biológické riziká/
Biolóskli rizici/Biyolojik riskler

ORIG HUM

Human/C човешки производ/Lidské/
Humanit/Human/δείγματα αναφοράς/
Humano/Inimpăritolu/Humaine/Ljudskog
porjekla/Humán/Origine Umana/
Žmogaus kilmės/Cilvēku izcelsmes/
Human/Menneske/Ludzka/Humano/
Origine umana/Человеческого
происхождения/Human/Ludské/
Humanega izvora/Ljudskog porekla/İnsan

ORIG MOU

From mouse/C миши производ/Мыši/
Fra mus/Maus/από ποντίκι/de ratón/
Hiirtelt/De souris/Mišijeg porjekla/
Egérbőli/Murino/Pelės kilmės/No peles/
Van muizen/Fra mus/Mysia/Do rato/De
la șoareci/Мышиного происхождения/
Från mus/Myšie/Mišjega izvora/Mišijeg
porekla/Farenden

ORIG BOV

Bovine/C говежди производ/Hovēži/
Bovin/Rind/από βοοειδή/Bovino/
Veistelt/Bovine/Rogate stoke/Szar-
vasmarha/Bovino/Jaučio/No liellopa/
Bovien/Bovin/Wolowy/Orwigo/Origine
bovină/крупного рогатого скота/Frán
ko/Hovädzie/Govejega izvora/Rogate
krupne stoke/Bovin

Reconstitute with/Разтваряне с/
Rozfedte pomocí/Rekonstitueres med/
Rekonstituieren mit/Ανασύσταση με/
Reconstituir con/Lahjendamine/
Reconstituer avec/Rekonstituiraite s/
Feloldáshoz/Ricostituire con/Atkurti,
ištirpdant su/Atšķaidīt ar/Reconstitue
met/Rekonstitueres med/Odtworzyć
za pomocą/Reconstituir com/A
se reconstitui cu/Разтворить в/
Rekonstituera med/Rozriedte pomocou/
Rekonstituiraite z/s/Ponovino formiranje
sa/Yeniden oluşturulur

Manufacturer/Производител/Výrobce/
Producent/Hersteller/Κατασκευαστής/
Fabricante/Tootja/Fabricant/Proizvođač/
Gyártó/Fabbricante/Gamintojas/
Ražotājs/Fabrikant/Produsent/
Producent/Fabricante/Producător/
Производитель/Tilverkare/Výrobca/
Izdelovalec/Proizvođač/Üretici

INSTRUCTIONS FOR USE

EN

INSTRUCTIONS FOR USE

Please visit our website www.fdi.com/ifu to obtain the Instructions For Use (IFU) in additional languages.

To ensure that you download the correct IFU for your kit lot, please select the revision corresponding to the issue date printed on the front page of the IFU provided with this kit.

Please follow the IFU carefully. Instructions for safe handling are found in the WARNINGS AND PRECAUTIONS section. Material Safety Data Sheets (MSDS) are available on our website www.fdi.com. If you do not have access to the internet, please contact your local distributor, or Fujirebio Diagnostics AB for assistance.

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BIBLIOGRAFIE

1. Bařt R., Feeney M., Lazarus H., Nadler L., Colvin R., Knapp R., (1981). Reactivity of a monoclonal antibody with human ovarian carcinoma. *J Clin Invest* 68: 1331 - 1337.
2. O'Brien T, Raymond L, Bannon G, Ford D, Hardardottir H, Miller F, Quirk J (1991). New monoclonal antibodies identify the glycoprotein carrying the CA125 epitope. *Am J Obstet Gynecol* 165: 1857 - 1864.
3. Nilsson O, Jansson E-L, Dahlen U, Nilsson K, Nuřtad K, Högberg T, Lindholm L (1994). In: Current Tumor Diagnosis: Applications, Clinical Relevance, Research, Trends. Ed R. Klapdor. Zucksverdt Verlag München, pp 401-405.
4. Nuřtad K. et al., (1996) Specificity and affinity of 26 monoclonal antibodies against the CA125 antigen: First report from the ISOBM TD-1 Workshop. *Tumor Biol* 17:196-219.
5. Bonfrer J., Duffy M., Radtke, M., Segurado O., Torre G., van Dalen A., Zwirner M. (1999) Tumor markers in gynaecological cancers-EGTM recommendations. *Anticancer Research* 19:2785-2820.
6. Tuxen M., Sölétormos, G., Ruřtin, G., Nelřtrop A., Dombernowsky, P. (2000) Biological variation and analytical imprecision of CA125 in patients with ovarian cancer. *Scand J Clin Invest* 60:713-722.
7. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
8. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).
9. Grimsey P et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *J Pharmacokinet.* 2017; 2(4); 247-256.
10. Jenkins Colon P, Greene D.N. Biotin Interference in Clinical Immunoassays. *JALM* 2018; 2(6); 941-951.
11. Chun KY. Biotin Interference in Diagnostic Tests. *Clinical Chemistry* 2017; 63; 619-620.

CanAg[®] is a registered trademark of Fujirebio Diagnostics AB

Fujirebio Diagnostics AB

Elof Lindälvs gata 13

SE-414 58 Göteborg

Sweden

Phone + 46 31 85 70 30

Fax + 46 31 85 70 40

info@fdab.com

www.fdab.com