Aesculap® Acculan 4

en Instructions for use/Technical description

Reciprocating saw GA334

USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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



Reciprocating saw GA334


Legend

- 1 Product (reciprocating saw)
- 2 Trigger (for oscillation frequency control)
- 3 Safety catch
- 4 Battery bay
- 5 Locking pin
- 6 Sterile funnel
- 7 Battery
- 8 Lid
- 9 Lid release element
- 10 Battery removal device
- 11 Saw blade
- 12 Rotating sleeve
- 13 Cut-outs (for centering sleeve)
- 14 Saw blade coupling
- 15 Cleaning template
- 16 Support foot
- 17 Rinsing adapter
- 18 Opening (at support foot)
- 19 Marking (at support foot)
- 20 Locking sleeve (to couple the support foot)
- 21 Centering pins
- 22 Inspection hole (for the bottom of the slot of the saw blade coupling)

The depictions are schematic only.

Symbols on product and packages

	<p>Caution Observe important safety information such as warnings and precautions in the instructions for use.</p>
 YYYY-MM	<p>Maintenance label Indication of the next maintenance appointment (Date: Year-Month)</p>
	<p>Two-dimensional machine readable code The code contains a unique serial number that can be used for electronically tracking single instruments. The serial number is based on the worldwide standard sGTIN (GS1).</p>
	<p>Manufacturer</p>
	<p>Date of manufacture</p>
	<p>Follow the instructions for use</p>
 	<p>Labeling of electrical and electronic devices pursuant to directive 2002/96/EG (WEEE)</p>

	<p>Classification Type BF</p>
	<p>Manufacturer's batch designation</p>
	<p>Manufacturer serial number</p>
	<p>Manufacturer order number</p>
	<p>Temperature limits during transport and storage</p>
	<p>Air humidity limits during transport and storage</p>
	<p>Atmospheric pressure limits during transport and storage</p>
S9	<p>Operating mode</p>
Rx only	<p>According to US federal law, this product may only be sold by a physician or by the order of a physician</p>
	<p>CE mark according to Regulation (EU) 2017/745</p>
	<p>Medical device</p>
	<p>Oscillation frequency control</p>
	<p>Locking sleeve open/closed</p>

Contents

1.	About this document	3
1.1	Scope	3
1.2	Warnings	3
2.	General information	3
2.1	Intended use	3
2.2	Main functions and design characteristics	3
2.2.1	Operating mode	3
2.3	Indications	4
2.4	Absolute contraindications	4
2.5	Relative contraindications	4
3.	Safe handling	4
4.	Product description	4
4.1	Scope of supply	4
4.2	Components required for operation	4
4.3	Operating principle	4
5.	Preparation	4
6.	Working with the device	4
6.1	System set-up	5
6.1.1	Connecting the accessories	5
6.1.2	Inserting the rechargeable battery	5
6.1.3	Intraoperative battery change	5
6.1.4	Removing the rechargeable battery	5
6.1.5	Protection against inadvertent activation	5
6.1.6	Attaching and detaching the saw blade	6
6.2	Function checks	6
6.3	Safe operation	6
6.3.1	Operating the product	6
7.	Validated reprocessing procedure	7
7.1	General safety notes	7
7.2	General information	7
7.3	Reusable products	7
7.4	Preparations at the place of use	7
7.5	Preparation before cleaning	7
7.6	Product-specific safety instructions for the processing procedure	7
7.7	Manual cleaning with wipe disinfection	8
7.8	Automatic cleaning/disinfection with manual pre-cleaning	9
7.8.1	Manual pre-cleaning with a brush	9
7.8.2	Mechanical alkaline cleaning and thermal disinfection	9
7.9	Inspection, maintenance and checks	10
7.10	Packaging	10
7.11	Steam sterilization	10
7.12	Storage	10
8.	Maintenance	10
9.	Troubleshooting list	11
10.	Technical Service	12
11.	Accessories/Spare parts	13
12.	Technical data	13
12.1	Classification according to Regulation (EU) 2017/745	13
12.2	Performance data, information about standards	13
12.3	Nominal operating mode	13
12.4	Ambient conditions	13
13.	Disposal	13

1. About this document

1.1 Scope

These instructions for use apply to the products listed below:

Art. no.	Name
GA334	Reciprocating saw

► For instructions for use of specific items and information on material compatibility, see B. Braun eIFU at eifu.bbraun.com

1.2 Warnings

Warnings indicate risks to patients, users, and/or the product that could develop while using the product. Warnings are marked as follows:

⚠ DANGER

Indicates a potential danger. If it is not prevented, death or severe injuries may result.

⚠ WARNING

Indicates a potential danger. Minor or moderate injuries may result if this is not prevented.

⚠ CAUTION

Indicates a potential damage to property. The product may become damaged if it is not prevented.

2. General information

2.1 Intended use

Task/Function

When combined with the appropriate tool, and where required with the corresponding protection guide, the Reciprocating saw GA334 is used to work on hard tissue, cartilage and related materials, as well as bone replacement materials.

Application Environment

The product fulfills the requirements for type BF pursuant to IEC/DIN EN 60601-1 and is used in operating rooms in sterile environments of explosion risk areas (such as areas with pure oxygen or anaesthesia gases).

2.2 Main functions and design characteristics

Oscillation frequency	min. 0 rpm to max. 15 000 rpm
-----------------------	-------------------------------

2.2.1 Operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

■ 30 second application, 30 second pause

■ 5 repetitions

■ 30 min cooling time

■ Max. Temperature 48 °C

Electrical systems generally heat up during continual operation. It is advised to give the system a break after use to cool down, as listed in the table on operating mode.

Heating depends on the tool used and the load. After a certain number of repetitions, the system should cool down. This procedure prevents the system overheating as well as possible injury to the patient or user.

The user is responsible for the use and adherence to the pause sequence described.

2.3 Indications

Type and area of application depend on the selected tool.

2.4 Absolute contraindications

The product is not approved for use in the central nervous system or central circulatory system.

2.5 Relative contraindications

The safe and effective use of the product depends greatly on factors that only operators themselves can control. The information provided therefore only represents general conditions.

The clinically successful use of the product depends on the knowledge and experience of the surgeon. They must decide which structures can be sensibly treated and must take the safety information and warnings in the instruction for use into account.

3. Safe handling

WARNING

Risk of injury and material damage if this product is not used as intended!

- ▶ Use the product only for its intended purpose.

WARNING

Risk of injury and damage to property due to improper handling of the product!

- ▶ Follow the instructions for use of all products used.
- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Check the functionality and proper condition before using the product.
- ▶ Observe "Guidance on Electromagnetic Compatibility (EMC) for Acculan Components" TA022450, see B. Braun eIFU at eifu.bbraun.com
- ▶ To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Keep the instructions for use accessible for the user.
- ▶ Comply with valid standards.
- ▶ Ensure that the electrical installation of the room is consistent with the requirements of IEC/DIN EN.
- ▶ Do not operate the product in explosion-hazard areas.
- ▶ Sterilize product before use.
- ▶ When using ECCOS holder systems, observe relevant instructions for use TA009721, see B. Braun eIFU at eifu.bbraun.com

Note

The user must report all serious incidents occurring in connection with the product to the manufacturer and to the competent authority of the country in which the user is registered.

4. Product description

4.1 Scope of supply

Art. no.	Designation
GA334	Reciprocating saw
GA674890	Rinsing adapter installed
GA678	Sterile funnel
GB440215	Cleaning template
TA014545	Instructions for use of reciprocating saw GA334 (flyer)

4.2 Components required for operation

- NiMH battery long GA676 (charged)
- Sterile funnel GA678
- Lid GA675
- Support foot GB436R (as required)
- Saw blade

4.3 Operating principle

The product **1** contains an electric motor, which is powered by a replaceable battery **7**.

The charged non-sterile battery **7** is inserted into the product **1** via the sterile funnel **6** and sealed with the lid **8** so that it is sterile.

The oscillation frequency is controlled electronically and can be continuously regulated with the trigger **2**.

The product **1** has a coupling **14**, which allows for the attachment of various saw blades **11**. These saw blades **11** lock automatically when attached to the product **1**.

The saw blade **11** can be removed by activating a rotating sleeve **12**.

The product **1** has an interface, which allows for the attachment of support foot **16**.

5. Preparation

Aesculap assumes no liability if the following rules are not followed:

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Prior to use, inspect the product and its accessories for any visible damage.
- ▶ Use the products and their accessories only if they are in perfect technical condition.

6. Working with the device

WARNING

Risk of infection and contamination!

Product is delivered unsterilized!

- ▶ Sterilize the product before use pursuant to the operating instructions.

WARNING

Risk of injury and material damage due to accidental activation of the product!

- ▶ Products which are not being actively used must be secured against accidental activation (position OFF).

⚠ WARNING

Risk of injury and material damage due to inappropriate use of tools!

- ▶ Always follow the safety advice and information given in the instructions for use.
- ▶ When coupling/uncoupling, handle tools with cutting edges with care.

⚠ WARNING

Damage to the product if dropped!

- ▶ Use the products only if they are in perfect technical condition, see Function check

⚠ WARNING

Risk of burns to skin and tissue through blunt tools or if product has not been maintained properly!

- ▶ Use tools only if they are in perfect condition.
- ▶ Replace blunt tools.
- ▶ Maintain the product properly, see maintenance guide.

6.1 System set-up**6.1.1 Connecting the accessories**

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

All configurations must fulfill the fundamental standard IEC/DIN EN 60601-1. The person connecting the devices with each other is responsible for the configuration and must ensure that the fundamental standard IEC/DIN EN 60601-1 or relevant national standards.

- ▶ Follow the instructions for use of individual accessories.
- ▶ In case of questions, contact your B. Braun/Aesculap partner or Aesculap Technical Service, address see Technical Service.

6.1.2 Inserting the rechargeable battery**⚠ CAUTION**

The product does not function when using the wrong battery!

- ▶ Do not use battery GA666 (with red base) or GA346.
- ▶ Product 1 with the battery shaft 4 must be turned upwards and the sterile funnel 6 (sterile) attached, see Fig. A.
- ▶ Battery 7 (non-sterile) is inserted into the battery shaft 4 by a second (non-sterile) person, see Fig. A.

Note

After the battery has been inserted, multiple signals will sound, which indicate that the product is ready to use.

- ▶ After insertion of the battery, have the sterile funnel 6 (non-sterile) removed by a second person.
- ▶ Attach the lid 8 (sterile) in such a manner that it clicks into place with both release catches 9.

Note

The sterility of the product is only guaranteed with a correctly placed lid.

6.1.3 Intraoperative battery change**⚠ WARNING**

Risk of injury from hot battery!

Battery may be hot after use in the machine.

- ▶ Remove the battery with the battery removal device and let it cool.

The battery removal device is used for replacing the battery during surgery while ensuring that sterile conditions are maintained.

- ▶ Turn the product 1 with the battery shaft 4 upwards.
- ▶ Press both safety catches 9 on the lid 8 at the same time and remove the lid 8.
- ▶ Attach the sterile battery removal device 10, see Fig. B.
- ▶ Gently shake the product 1 with the attached battery removal device 10 with battery compartment 4 facing downwards.
The battery 7 slides easily into the battery removal device 10.
- ▶ Hand battery removal device 10 including dead battery 7 to a non-sterile person.
- ▶ Insert a charged battery 7, see Inserting the rechargeable battery.

6.1.4 Removing the rechargeable battery**⚠ WARNING**

Risk of injury from hot battery!

Battery may be hot after use in the machine.

- ▶ Let the battery cool in the machine and then remove.
– or –
- ▶ Remove the battery with the battery removal device and let it cool.

⚠ CAUTION

Damage to the battery caused by knocking it against hard objects!

- ▶ Only remove battery by tapping the product on the palm of the hand.

⚠ CAUTION

Damage to, or destruction of the batteries due to processing!

- ▶ Do not sterilize batteries.

Rechargeable battery must be removed after each surgical procedure and prior to processing of the device.

Note

The battery removal device can be used for easier battery removal, see Fig. B.

- ▶ Turn the product 1 with the battery shaft 4 upwards.
- ▶ Press both safety catches 9 on the lid 8 at the same time and remove the lid 8.
- ▶ Tightly grasp the battery shaft 4 at the lower end.
- ▶ Tap the battery shaft 4 on the palm of the hand until the battery 7 slides out of the battery shaft 4 and can be removed.

6.1.5 Protection against inadvertent activation

To prevent inadvertent activation of the product during tool change, the trigger can be locked.

To lock trigger:

- ▶ Twist the safety catch 3 to position OFF.
The trigger 2 is blocked and the product 1 cannot be operated.

To unlock trigger:

- ▶ Twist the safety catch 3 to position ON.
The trigger 2 is unlocked and the product 1 can be operated.

6.1.6 Attaching and detaching the saw blade

⚠ WARNING

Risk of injury when attaching/removing sawblades or support foot in the ON position through inadvertent activation of the product!

- ▶ Only attach/remove saw blades or support foot in the OFF position.

Coupling the saw blade

- ▶ Secure product 1 against inadvertent activation with the safety catch 3.
- ▶ Turn rotating sleeve 12 as far as it will go against the direction indicated by the arrow.
- ▶ Insert saw blade 11 for pulling or pushing incision in saw blade coupling 14.
Saw blade 11 engages.

Uncoupling the saw blade

- ▶ Secure product 1 against inadvertent activation with the safety catch 3.
- ▶ Turn rotating sleeve 12 180° as far as it will go in the direction indicated by the arrow.
- ▶ Pull out saw blade 11 from the saw blade coupling 14.

Coupling the support foot

- ▶ Secure product 1 against inadvertent activation with the safety catch 3.
- ▶ Open locking sleeve 20 at support foot 16. Marking (on support foot) 19 points to symbol "open."
- ▶ Pull back and hold locking sleeve 20 towards support foot 16.
- ▶ Push the support foot 16 all the way over the rotating sleeve 12 of the saw blade coupling 14 so that the saw blade 11 grips into the opening (on the support foot) 18 and the centering pins 21 sit in the cut-outs 13. The teeth of saw blade 11 must point away from support foot 16.
- ▶ Release and turn locking sleeve 20 until marking (at support foot) 19 points to symbol "closed."
The support foot 16 is locked in place.

Uncoupling the support foot

- ▶ Secure product 1 against inadvertent activation with the safety catch 3.
- ▶ Open locking sleeve 20 at support foot 16.
The support foot 16 is unlocked.
- ▶ Push the support foot 16 in the direction of the product 1 and pull the locking sleeve 20 from the product 1.

6.2 Function checks

The function checks must be carried out prior to each use and after each intraoperative battery change.

- ▶ Make certain there is a battery in the battery compartment.
- ▶ Ensure that the lid has completely clicked into place.
- ▶ Pull on the saw blade to check that it is securely coupled.
- ▶ Check that the cutting edges of the saw blade do not exhibit any mechanical damage.
- ▶ Check secure seating of the support foot: ensure locking sleeve is closed.
- ▶ Release product for use (ON position).
- ▶ Briefly operate product at maximum oscillation frequency.
- ▶ Do not use the product if it is damaged or defective.
- ▶ Set aside the product if it is damaged.

6.3 Safe operation

⚠ WARNING

Coagulation of patient tissue or risk of burns for patients and user through hot product!

- ▶ Cool the tool during operation.
- ▶ Keep product/tool out of the reach of patients.
- ▶ Let the product/tool cool down.
- ▶ Use a cloth to protect against burns when changing the tool.

⚠ WARNING

Risk of infection from aerosol formation!

Risk of injuries caused by particles coming loose from the tool!

- ▶ Use suitable protection (such as waterproof protective clothing, face mask, safety gases, suction).

⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Always carry out a function check prior to using the product.

⚠ WARNING

Risk of injury when using the product beyond the field of view!

- ▶ Apply the product only under visual control.

⚠ WARNING

Risk of injury and damage to the tool/system!

The tool may get caught in drapes (such as textiles).

- ▶ Do not let the tool come into contact with drapes (such as textiles) during operation.

Note

The product drive motor is operated using a magnetic sensor system. In order to prevent inadvertent activation of the motor, the product may not be exposed to any magnetic fields (such as magnetic instrument pads).

6.3.1 Operating the product

Note

The slight whistling noise when the product is starting up is due to the construction.

- ▶ Activate trigger 2.
The oscillation frequency of the product 1 is continuously variable.

7. Validated reprocessing procedure

7.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

See also B.Braun eIFU at eifu.bbraun.com for current information on processing and compatibility of material

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

On stainless steel, residues containing chlorine or chloride (such as surgical residues, drugs, saline solutions in water for cleaning, disinfection and sterilization) may lead to corrosion (pitting corrosion, tensile corrosion) and thus to the destruction of the product. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Material damage (such as corrosion, cracks, breaks, premature aging or swelling).
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org link to "AKI-Brochures", "Red brochure".

7.3 Reusable products

The service life of the product is limited by damage, normal wear, type and duration of the application, handling, storage and transportation of the product.

Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product.

7.4 Preparations at the place of use

- ▶ Remove all attached components from the product (tool and accessories).
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

7.5 Preparation before cleaning

- ▶ Prior to first mechanical cleaning/disinfection: Assemble ECCOS holders in a suitable basket (e.g. GB243800).
- ▶ Insert product 1 in the correct position into the ECCOS holders, see Fig. B.

7.6 Product-specific safety instructions for the processing procedure

⚠ CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for plastic material and high-grade steel,
 - and which do not affect softeners (e.g., in silicone).
- ▶ Do not use cleaning agents that contain acetone.
- ▶ Observe information concerning concentration, temperature and exposure time.
- ▶ Do not exceed the maximum temperature of 60 °C in chemical cleaning and/or disinfection.
- ▶ Do not exceed maximum temperature during thermal disinfection using 96 °C DI water.
- ▶ Dry the product for at least 10 minutes at a maximum of 120 °C.

⚠ CAUTION

Damage or destruction of the batteries due to processing!

- ▶ Protect battery from moisture.

Note

The drying time listed is only indicative. It has to be checked taking into account the specific conditions (e.g. load) and adjusted if necessary.

7.7 Manual cleaning with wipe disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Pre-cleaning	RT (cold)	≥2	–	D–W	Until visually clean
II	Cleaning with enzyme solution	RT (cold)	≥2	0.8	D–W	pH neutral*
III	Intermediate rinse	RT	≥5	–	D–W	–
IV	Drying	RT	–	–	–	–
V	Wipe disinfection	–	>1	–	–	Meliseptol HBV wipes with 50 % propan-1-ol
VI	Final rinse	RT (cold)	0.5	–	FD–W	–
VII	Drying	RT	–	–	–	–

D–W: Drinking water

FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

* Suitable enzyme solution: Helizyme, Cidezyme (the latter used for validation purposes)

- ▶ Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.
- ▶ Push the saw blade coupling **14** towards the housing as far as it will go.

Phase I

- ▶ Mobilize non-rigid components during cleaning, e.g. pull locking sleeve **20** back from the support foot **16**.
- ▶ Clean the product under running water, using a suitable plastic cleaning brush until all visible residues have been removed from the surfaces.
- ▶ Use cleaning template **15** to clean the slot of the saw blade coupling **14**. Insert cleaning template into the slot and carefully push out the debris through the slot or the inspection holes at least 5 times in both directions. Turn template 180° in order to do so.
- ▶ Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.
- ▶ Use rinsing adapter **21** to clean the inside of the saw blade coupling **14**:
 - Fit rinsing adapter **21** onto rotating sleeve **12**.
- ▶ Using an appropriate syringe, rinse slot manually with 100 ml tap water.
- ▶ Remove rinsing adapter.

Note

For details on difficult to access surfaces, see *Acculan pre-cleaning and care information TA016000* (available in Aesculap Extranet at <https://extranet.bbraun.com>).

Phase II

- ▶ Pull back locking sleeve **20** on support foot **16**, push cleaning template **15** between sleeve and support foot and release locking sleeve **20**.
- ▶ Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- ▶ Spray products with a pH neutral enzyme solution, let soak in for at least 2 minutes and then wipe off.
- ▶ Contamination should be removed with a lint-free cloth or soft brush moistened with enzyme cleaner.
- ▶ Flexible components (such as sleeves) should be rinsed for 20 seconds with the water pistol (cold water, at least 2.5 bar).
- ▶ After manual cleaning, check visible surfaces and areas of flexible components for residues.

Phase III

- ▶ Rinse product under running tap water for at least 5 minutes.
- ▶ Flexible components (such as sleeves) will move during cleaning.
- ▶ If necessary, repeat the cleaning process (phase 1 to 3).
- ▶ Remove cleaning template **15** from support foot **16**.

Phase IV

- ▶ Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

Phase V

- ▶ Wipe all surfaces of the product with a single-use disinfecting wipe.

Phase VI

- ▶ Rinse disinfected surfaces after the prescribed reaction time for at least 1 minute under running demineralized water.
- ▶ Drain any remaining water fully.

Stage VII

- ▶ Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

7.8 Automatic cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must fundamentally have a tested efficacy (such as FDA approval or CE label pursuant to DIN EN ISO 15883).

Note

The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

7.8.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemistry/Note
I	Rinse	RT (cold)	–	–	D–W	Until visually clean
II	Brush	RT (cold)	–	–	D–W	until visually clean

D–W: Drinking water
RT: Room temperature

- ▶ Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.
- ▶ Push the saw blade coupling **14** towards the housing as far as it will go.

Phase I

- ▶ Mobilize non-rigid components during cleaning, e.g. pull locking sleeve **20** back from the support foot **16**.
- ▶ Thoroughly clean the product under running water.
- ▶ Use cleaning template **15** to clean the slot of the saw blade coupling **14**. Insert cleaning template into the slot and carefully push out the debris through the slot or the inspection holes at least 5 times in both directions. Turn template 180° in order to do so.
- ▶ Manually clean the contact sites of the rinsing adapter GA674890 and the rotating sleeve **12** thoroughly.

Phase II

- ▶ Mobilize non-rigid components during cleaning, e.g. pull locking sleeve **20** back from the support foot **16**.
- ▶ Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.
- ▶ After manual preparation, check visible surfaces for residue and repeat the pre-cleaning process as needed.

Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available in Aesculap Extranet at <https://extranet.bbraun.com>).

7.8.2 Mechanical alkaline cleaning and thermal disinfection

Type of device: Single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	<ul style="list-style-type: none"> ■ Concentrate, alkaline: <ul style="list-style-type: none"> – pH ~ 13 – <5 % anionic surfactant ■ Working solution 0.5% <ul style="list-style-type: none"> – pH ~ 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfection	90/194	5	FD–W	–
V	Drying	max. 120/248	min. 10 min	–	–

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Insert the product in its correct position into the ECCOS holder.
- ▶ Pull back locking sleeve **20** on support foot **16**, push cleaning template **15** between sleeve and support foot and release locking sleeve **20**.
- ▶ Then connect the rinsing adapter GA674890 to the rotating sleeve **12** and connect to the rinse connector of the cleaning/disinfection device/rinse cart.
- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.
- ▶ Remove cleaning template **15** from support foot **16**.
- ▶ Check the floor of the saw blade coupling **14** slit (open inspection holes) for cleanness.
- ▶ Repeat the cleaning/disinfection process if necessary.

7.9 Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.

Note

Aesculap additionally recommends the occasional spraying of moving parts (such as knobs, connections) with STERILIT Power Systems oil spray.

- ▶ Inspect the product after every cleaning and disinfecting cycle for: cleanliness, damage, function, abnormal operation noise, excessive heat or strong vibration.
- ▶ Inspect the saw blades for broken, damaged or blunt cutting edges.
- ▶ Set aside the product if it is damaged.

7.10 Packaging

- ▶ Observe the instructions for use of the respective packaging and holders (e.g. Instructions for use TA009721 for Aesculap ECCOS holder system).
- ▶ Insert product **1** in the correct position in the ECCOS holders, see Fig. B.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

7.11 Steam sterilization

⚠ CAUTION

Damage or destruction of the batteries due to processing!

- ▶ **Do not sterilize the rechargeable battery.**

Note

Remove all attached components from the product (tools, accessories) before sterilization.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Use a validated sterilization method:
 - Steam sterilization in fractionated vacuum process
 - Steam sterilizer in accordance with DIN EN 285 and validated in accordance with DIN EN ISO 17665
 - Sterilization in fractionated vacuum process at 134°C, holding time 5 min

When sterilizing multiple products in one steam sterilizer:

- ▶ Ensure that the maximum permissible load of the steam sterilizer according to the manufacturer's specifications is not exceeded.

7.12 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

8. Maintenance

To ensure reliable operation, the product must be maintained in accordance with the maintenance labeling or at least once a year.

YYYY-MM

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

9. Troubleshooting list

► Have defective products repaired by Aesculap Technical Service, see Technical Service.

Fault	Cause	Detection	Remedy
Product not running	No battery	No battery in the battery shaft	Insert battery.
	Battery not charged	No signal tone when inserting the battery	Charge battery in charger.
	Battery broken	No signal tone when inserting the battery	Have the manufacturer repair the battery.
	Product in secured OFF position	Safety catch is in OFF position	Turn the safety catch into the ON position.
	Product defective	Product not running	Have the manufacturer repair the product.
Product becomes too hot	Excessive use	Product heating	Observe instructions for use (nominal operating mode).
	Reprocessing/maintenance carried out incorrectly	Product heating	Observe instructions for use (processing, care). Preventive measure: Oil product before every sterilization.
	Fall damage, product defective	Product heating	Have the manufacturer repair the product.
	Blunt tool	Tool and product become hot	Change tool.
Insufficient power	Product broken	Insufficient product performance	Observe instructions for use (processing, care). Preventive measure: Oil product before every sterilization. Have the manufacturer repair the product.
		Steep rise in temperature after a short period of time	Observe instructions for use (nominal operating mode). Have the manufacturer repair the product
	Blunt tool	Tool blades worn	Change tool.
Loud running noise	Broken gearing/ball bearings of the product	Loud, perceptible noise during operation	Follow operating instructions (preparation, care). Preventative: lubricate product before every sterilization. Have the manufacturer repair the product
Lid cannot be attached/removed	Lid not compatible	Lid does not click into place	Use the compatible lid for GA334.
	Lid deformed/broken	Lid is difficult or cannot be mounted/removed	Have the manufacturer repair the lid.
	Connection on product deformed/broken	Lid is hard or impossible to install/remove	Have the manufacturer repair the product.
Battery cannot be inserted/removed	Battery not compatible	Battery cannot be completely inserted into the battery compartment	Use the compatible battery for GA334.
	Battery deformed/broken	Battery is hard or impossible to insert/remove	Have the manufacturer repair the battery.
	Battery shaft to product deformed/defective	Battery is difficult or cannot be mounted/removed	Have the manufacturer repair the product
Sterile funnel or battery removal device cannot be attached	Sterile funnel or battery removal device not compatible	Sterile funnel or battery removal device not attachable to the battery shaft	Use suitable sterile funnel or battery removal device for GA334.
	Sterile funnel or battery removal device deformed/broken	Sterile funnel or battery removal device difficult or not attachable	Change sterile funnel or battery device.
	Battery compartment at product deformed/broken	Sterile funnel or battery removal device is hard or impossible to attach	Have the manufacturer repair the product.

Fault	Cause	Detection	Remedy
Speed control knob cannot be actuated	Product in locked OFF position	Safety catch is in OFF position	Turn the safety catch to the ON position.
	Speed control knob sticks/is defective	Speed control knob cannot be actuated	Have the manufacturer repair the product.
Saw blade does not move	Noisy motor	Broken gearings	Have product repaired by the manufacturer
Insufficient cutting performance of the saw blade	Saw blade teeth worn	Blunt saw blade	Replace saw blade
	Battery power is low	Power/speed of drive inadequate	Charge battery
	Battery empty/broken Symbol "Battery replacement recommended" illuminated on the charger		Replace battery
The saw blade cannot be coupled	Debris on the saw blade or in the coupling	Debris	Clean saw blade/coupling with cleaning template 7
	Deformed coupling parts/saw blades	Deformation	Have product repaired by the manufacturer
	-	Incompatible saw blade	Insert a compatible Aesculap saw blade, see Accessories/Spare parts
Saw blade cannot be locked and can be pulled out during function check	Saw blade can be pulled out without rotating sleeve action	Rotating sleeve in incorrect position (not returning automatically to rest position)	Clean clutch if soiled If damaged, have device repaired by the manufacturer
Saw blade cannot be uncoupled	-	Rotating sleeve in incorrect position	Turn rotating sleeve 180° in direction of the arrow, see Attaching and detaching the saw blade
	Saw blade coupling warped or bent saw blade	Mechanical damage	Have repaired by manufacturer
Support foot cannot be coupled	Support foot cannot be coupled	Saw blade bent	Replace saw blade
	Support foot cannot be coupled	Support foot bent	Have support foot repaired by the manufacturer
	Marking on support foot points to "Closed" symbol.	Locking sleeve is closed	Open locking sleeve, see Coupling the support foot.

10. Technical Service

DANGER

Danger to life of patients and users in case of malfunctions and/or failure of protective measures!

- ▶ Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.

CAUTION

Modifications of medical devices may result in a loss of potential guarantee/warranty claims and forfeiture of applicable licenses.

- ▶ Do not modify the product.
- ▶ Contact national B. Braun/Aesculap representative for service and repair.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95 -1601

Fax: +49 (7461) 14 -939

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

11. Accessories/Spare parts

Art. no.	Name
GA674890	Rinsing adapter
GA675	Lid
GA676	NiMH battery long
GA678	Sterile funnel
GA679	Battery removal device
GB243R GB244R	ECCOS basket with holder for Acculan 4
GB600	STERILIT Power Systems oil spray
GB436R	Support foot
GB440215	Cleaning plate
GB495R	ECCOS holder for machine
GB487R	ECCOS holder for lid
GB488R	ECCOS holder for sterile funnel
GB489R	ECCOS holder for battery removal device
GB686R	ECCOS support foot holder
TA014544	Instructions for use of reciprocating saw GA334 (A4 for ring binder)
TA014545	Instructions for use of reciprocating saw GA334 (flyer)

12. Technical data

12.1 Classification according to Regulation (EU) 2017/745

Art. no.	Designation	Class
GA334	Reciprocating saw	Ila

12.2 Performance data, information about standards

Max. power	approx. 250 W
Max. oscillation frequency	15 000 rpm
Weight (operational)	1.54 Kg ±10 %
Dimensions (L x W x H, ready for operation)	207 mm x 180 mm x 55 mm ±5 %
Application part	Type BF
EMC	IEC/DIN EN 60601-1-2
Conforming to standard	IEC/DIN EN 60601-1

After 500 processing cycles, the manufacturer tested the product and it passed.

12.3 Nominal operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

- 30 second application, 30 second pause
- 5 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

12.4 Ambient conditions

	Operation	Transport and storage
Temperature	10 °C to 27 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

13. Disposal

⚠ WARNING

Risk of infection from contaminated products!

- Observe national regulations when disposing of or recycling the product, its components, and their packagings.

Note

The operator must process the product before disposal, see Validated reprocessing procedure.

The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.) Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- For questions regarding the disposal of the product, please contact your local B. Braun/Aesculap representative, see Technical Service.

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




without base perforation = JK

with base perforation = JN

Lids for individual color combination

Inner lid

Outer lid

		Outside dimensions Length x Width x Height (mm)
JK339 without carrying handle		300 x 274 x 57
JK340	JN340	300 x 274 x 90
JK341	JN341	300 x 274 x 120
JK342	JN342	300 x 274 x 135
JK344	JN344	300 x 274 x 187
JK346	JN346	300 x 274 x 247

for JK339	for JK340-346	Color	for JK339-346	Color
JK395	JK385	red	JK375A	red
JK396	JK386	blue	JK376A	blue
JK397	JK387	green	JK377A	green
JK398	JK388	gold	JK378A	gold
JK399	JK389	silver	JK379A	silver

Dimensions complete - Length x Width x Height (mm)

	Outside dimensions Inside dimensions	Outside dimensions Inside dimensions
Body 1/2	with inner lid	with outer lid
JK339 without carrying handle	300 x 285 x 67 260 x 265 x 35	300 x 291 x 75 260 x 265 x 35
JK340	300 x 285 x 108 254 x 258 x 75	300 x 291 x 115 254 x 258 x 75
JK341	300 x 285 x 138 254 x 258 x 108	300 x 291 x 145 254 x 258 x 108
JK342	300 x 285 x 153 254 x 258 x 120	300 x 291 x 160 254 x 258 x 120
JK344	300 x 285 x 205 254 x 258 x 172	300 x 291 x 212 254 x 258 x 172
JK346	300 x 285 x 265 254 x 258 x 232	300 x 291 x 272 254 x 258 x 232

	Outside dimensions Inside dimensions	Outside dimensions Inside dimensions
Body 1/2	with inner lid	with outer lid
-		
JN340	300 x 285 x 108 254 x 258 x 63	300 x 291 x 115 254 x 258 x 63
JN341	300 x 285 x 138 254 x 258 x 93	300 x 291 x 145 254 x 258 x 93
JN342	300 x 285 x 153 254 x 258 x 108	300 x 291 x 160 254 x 258 x 108
JN344	300 x 285 x 205 254 x 258 x 160	300 x 291 x 212 254 x 258 x 160
JN346	300 x 285 x 265 254 x 258 x 220	300 x 291 x 272 254 x 258 x 220

All dimensions: Length x Width x Height (mm)


Fixation pins for ring handle instruments

	Diameter	height 29 mm
JF346	25 mm	

	Diameter	height 47 mm
JF343	6 mm	

Silicone mats (mesh)

	Outside dim.	slip-free storage of instruments good drying characteristics, gentle storage during mechanical processing. <i>Unsuitable for hot air sterilization!</i>
JF938	240 x 250 mm	
JF939	480 x 250 mm	
JF940	402 x 250 mm	
JF941	536 x 250 mm	
JF942	191 x 233 mm	
JF943	433 x 233 mm	
JF946	257 x 114 mm	

	Outside dim.	Positioning mat Specifications as above, but trimmed and perforated. <i>Unsuitable for hot air sterilization!</i>
JF934	276 x 125 x 17	

Silicone mats with perforation (nubs)

	Outside dim.	slip-free storage of instruments good drying characteristics, gentle storage during mechanical processing. Fitting to JF159R
JF944	248 x 102 mm	

	Outside dim.	Fitting to 1/2 basket
JF945	248 x 257 mm	

Unsuitable for hot air sterilization!

Reusable filters

JK090

Outside dim.

Ø 190

Aesculap reusable filter, round for use in up to 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 10

JK089

Outside dim.

225 x 225

Aesculap reusable filter, square for legacy Standard containers, service life 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 10

JK091

Outside dim.

95 x 215

Aesculap reusable filter for Mini and Dental containers **with filter holder JK098**, service life 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 2

JK092

Outside dim.

171 x 231

Aesculap reusable filter for Dental containers **without filter holder JK098**, service life 1,000 cleaning and sterilization cycles.

Order unit (PAK): Pack of 2

Marking pens

JK086

Marking pen for Aesculap reusable filters

Order unit (PAK): Pack of 3

All dimensions: Length x Width x Height (mm)