



en

AST2

04S90

G93275R03

B4S900

# Aspartate Aminotransferase2

FOR USE WITH

ARCHITECT

Revised July 2021.

REF 04S9020

REF 04S9030

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

For laboratory professional use only.

## NAME

Aspartate Aminotransferase2 (also referred to as AST2)

## INTENDED USE

The Aspartate Aminotransferase2 assay is used for the quantitation of aspartate aminotransferase (AST) in human serum or plasma on the ARCHITECT c Systems.

The Aspartate Aminotransferase2 assay is to be used as an aid in the diagnosis and treatment of certain liver diseases.

## SUMMARY AND EXPLANATION OF THE TEST

AST is a liver enzyme that is also found in the heart, skeletal muscle and kidney. AST has two isoenzymes that are found in either the mitochondria or cytoplasm in cells. Alcohol-induced hepatocyte injury induces predominantly mitochondrial damage. AST requires vitamin B6 as a cofactor for the enzymatic reaction. Decreased AST levels may indicate vitamin B6 deficiency and uremia.<sup>1-3</sup>

AST is most commonly used in conjunction with other laboratory findings [alanine aminotransferase (ALT) or lactate dehydrogenase (LDH)]. The DeRitis ratio (AST/ALT quotient) is found to be elevated in alcohol-induced liver disease, viral hepatitis, cirrhosis and acute fulminant hepatic failure. A ratio of greater than two suggests that alcohol is the cause of liver injury.<sup>1-3</sup>

## PRINCIPLES OF THE PROCEDURE

The Aspartate Aminotransferase2 assay is an automated clinical chemistry assay.

AST present in the sample catalyzes the transfer of the amino group from L-aspartate to  $\alpha$ -ketoglutarate, forming oxaloacetate and L-glutamate. Oxaloacetate in the presence of nicotinamide adenine dinucleotide (NADH) and malate dehydrogenase (MDH) is reduced to L-malate. In this reaction, NADH is oxidized to  $\text{NAD}^+$ . The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm due to the oxidation of NADH to  $\text{NAD}^+$ .

Methodology: NADH (without P-5'-P)

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

## REAGENTS

### Kit Contents

Aspartate Aminotransferase2 Reagent Kit 04S90

NOTE: Some kit sizes may not be available. Please contact your local distributor.

Volumes (mL) listed in the following table indicate the volume per cartridge.

REF	04S9020	04S9030
Tests per cartridge set	300	990
Number of cartridge sets per kit	4	4
Tests per kit	1200	3960
<b>R1</b>	28.2 mL	84.5 mL
<b>R2</b>	18.5 mL	53.8 mL
<b>R1</b> Active ingredients: L-aspartic acid (103.860 g/L), $\beta$ -NADH (0.610 g/L), lactate dehydrogenase (4.000 KU/L), and malate dehydrogenase (2.000 KU/L). Preservative: sodium azide.		
<b>R2</b> Active ingredient: $\alpha$ -ketoglutaric acid (6.570 g/L). Preservative: ProClin 300.		

## Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use
- **Rx ONLY**

## Safety Precautions

**CAUTION:** This product requires the handling of human specimens. It is recommended that all human-sourced materials and all consumables contaminated with potentially infectious materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate regional, national, and institutional biosafety practices should be used for materials that contain, are suspected of containing, or are contaminated with infectious agents.<sup>4-7</sup>

The following warnings and precautions apply to: <b>R1</b>	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

The following warnings and precautions apply to: <b>R2</b>	
<b>WARNING</b>	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
H402*	Harmful to aquatic life.
H412	Harmful to aquatic life with long lasting effects.
<b>Prevention</b>	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
P273	Avoid release to the environment.

Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

\* Not applicable where regulation EC 1272/2008 (CLP) has been implemented.

Follow local chemical disposal regulations based on your location along with recommendations and content in the Safety Data Sheet to determine the safe disposal of this product.

For the most current hazard information, see the product Safety Data Sheet.

Safety Data Sheets are available at [www.corelaboratory.abbott](http://www.corelaboratory.abbott) or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

### Reagent Handling

- Do not pool reagents within a kit or between kits.
- Do not use components from one lot with components from another lot.
- Do not reuse containers, caps or plugs due to the risk of contamination and the potential to compromise reagent performance.
- When either the **R1** or **R2** reagent cartridge becomes empty, replace both cartridges.
- Upon receipt, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

### Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
<b>Unopened</b>	2 to 8°C	Until expiration date	Store in upright position.
<b>Onboard</b>	System Temperature	30 days	
<b>Opened</b>	2 to 8°C	Until expiration date	Store in upright position.

Reagents may be stored on or off the ARCHITECT c System. If reagents are removed from the system, store at 2 to 8°C (with replacement caps) in their original boxes.

For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

### Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

## INSTRUMENT PROCEDURE

The Aspartate Aminotransferase2 assay file must be installed on the ARCHITECT c System prior to performing the assay.

Installation of all the required SmartWash updates on either the MULTIGENT Assay Disk Version 9.00 (or higher) or the Special Chemistry Assay Disk Version 7.00 (or higher) must be completed prior to performing the assay. See below for impacted assays:

Assay Name	Short Name	REF	Assay Number	VERSION	
				Conventional Units / Alternate Units	SI Units / Alternate Units
Benzodiazepines Qualitative	BenzQ	3L39	2853	5	5
Benzodiazepines Semiquantitative	BenzSQ	3L39	2852	5	5
Gentamicin	Gent	1E11	2867	9	9

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

### Alternate Result Units

Conversion formula:

(Concentration in Default result unit) x (Conversion factor) = (Concentration in Alternate result unit)

Default Result Unit	Conversion Factor	Alternate Result Unit
U/L	0.01667	μkat/L

## SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

### Specimen Types

The specimen types listed below were verified for use with this assay.

Other specimen types and collection tube types have not been verified with this assay.

Specimen Types	Collection Tubes
Serum	Serum Serum separator
Plasma	Lithium heparin Lithium heparin separator Sodium heparin

- Liquid anticoagulants may have a dilution effect resulting in lower concentration values for individual specimens.

The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

## Specimen Conditions

- Do not use:
  - heat-inactivated specimens
  - pooled specimens
  - grossly hemolyzed specimens
  - specimens with obvious microbial contamination
  - specimens with fungal growth
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

## Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Recentrifugation of Specimens

- Transfer specimens to a centrifuge tube and centrifuge.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

## Specimen Storage

Specimen Type	Temperature	Maximum Storage Time
Serum/Plasma	Room temperature (20 to 25°C)	4 days <sup>8, 9</sup>
	2 to 8°C	7 days <sup>8, 10</sup>
	-20°C	12 weeks <sup>11</sup>

Avoid multiple freeze/thaw cycles.<sup>11</sup>

It is the responsibility of the individual laboratory to determine specific specimen stability criteria for their laboratory per their laboratory workflow.

For additional information on sample handling and processing, refer to CLSI GP44-A4.<sup>12</sup> The storage information provided here is based on references.

Each laboratory may establish a range around -20°C from either the freezer manufacturer's specifications or your laboratory standard operating procedure(s) for specimen storage.

Stored specimens must be inspected for particulates. If present, mix with a low speed vortex or by inversion and centrifuge the specimen to remove particulates prior to testing.

## Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

Do not exceed the storage limitations listed above.

## PROCEDURE

### Materials Provided

04S90 Aspartate Aminotransferase2 Reagent Kit

### Materials Required but not Provided

- Aspartate Aminotransferase2 assay file found on [www.corelaboratory.abbott](http://www.corelaboratory.abbott)
- 04V1501 Consolidated Chemistry Calibrator, if using the Calibration method
- Controls containing aspartate aminotransferase
- Saline (0.85% to 0.90% NaCl) for specimen dilution

For information on materials required for operation of the instrument, refer to the ARCHITECT System Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

### Assay Procedure

For a detailed description of how to run an assay, refer to the ARCHITECT System Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the ARCHITECT System Operations Manual, Section 5 to ensure sufficient specimen is present.
- Minimum sample cup volume is calculated by the system and printed on the Order List report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Minimum sample volume requirements:
  - Sample volume for single test: 5.3 µL.  
NOTE: This amount does not include the dead volume plus the additional over-aspiration volume. For total sample volume requirements, refer to the ARCHITECT System Operations Manual, Section 5.
- Refer to the Consolidated Chemistry Calibrator package insert **REF** 04V1501 and/or commercially available control material package insert for preparation and usage.
- For general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

### Sample Dilution Procedures

Sample dilutions have not been evaluated for the Aspartate Aminotransferase2 assay. Samples with an aspartate aminotransferase value exceeding 4001 U/L (66.70 µkat/L) are flagged with code "> 4001 U/L" ("> 66.70 µkat/L").

For details on configuring automated dilutions, refer to the ARCHITECT System Operations Manual, Section 2.

### Calibration

For instructions on performing a calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Calibration can be performed using one of 2 methods:

- Calibration method, using the Consolidated Chemistry Calibrator **REF** 04V1501. For the Calibration method, use assay file AST2.
- Calibration Factor method, using a fixed calibration factor value to calculate the result. For the Calibration Factor method, use assay file AST2F.

Calibration is stable for approximately 30 days (720 hours) but is required with each change in reagent lot. Verify calibration with at least 2 levels of controls according to the established quality control requirements for your laboratory. If control results fall outside acceptable ranges, recalibration may be necessary.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

### Quality Control Procedures

As appropriate, refer to your laboratory standard operating procedure(s) and/or quality assurance plan for additional quality control requirements and potential corrective actions.

- At least two levels of controls (low and high) are to be run every 24 hours.
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

### Quality Control Guidance

Refer to "Basic QC Practices" by James O. Westgard, Ph.D. for guidance on laboratory quality control practices.<sup>13</sup>

## RESULTS

### Calculation

#### Calibration method

The Aspartate Aminotransferase2 (AST2) assay utilizes the Linear data reduction method to generate a calibration and results.

#### Calibration Factor method

The Aspartate Aminotransferase2 (AST2F) assay utilizes the Factor data reduction method to generate a calibration and results.

The calibration factors for the Aspartate Aminotransferase2 assay are 9246 (ARCHITECT c8000) and 8757 (ARCHITECT c4000 and ARCHITECT c16000).

The Aspartate Aminotransferase2 assay is traceable to the IFCC (International Federation of Clinical Chemistry) reference method.<sup>14</sup> The assigned values for the calibrator and the calibration factor are traceable to the standardization.

For additional information on results calculations, refer to the ARCHITECT System Operations Manual, Appendix C.

### Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

### Reportable Interval

Based on representative data for the limit of quantitation (LoQ) and the limit of detection (LoD), the ranges over which results can be reported are provided below according to the definitions from CLSI EP34, 1st ed.<sup>15</sup>

	U/L	μkat/L
Analytical Measuring Interval (AMI) <sup>a</sup>	5-4001	0.08-66.70
Reportable Interval <sup>b</sup>	4-4001	0.07-66.70

<sup>a</sup> AMI: The AMI extends from the LoQ to the upper limit of quantitation (ULoQ). This is determined by the range of values in U/L (μkat/L) that demonstrated acceptable performance for linearity, imprecision, and bias.

<sup>b</sup> The reportable interval extends from the LoD to the upper limit of the AMI.

NOTE: The default Low Linearity value of the assay file corresponds to the lower limit of the reportable interval of 4 U/L (0.07 μkat/L). To flag values using the lower limit of the analytical measuring interval of 5 U/L (0.08 μkat/L), the operator must edit the Low Linearity value.

For detailed information on editing the result settings of assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

### LIMITATIONS OF THE PROCEDURE

- Results should be used in conjunction with other data; e.g., symptoms, results of other tests, and clinical impressions.
- Substances that demonstrated interference with the Aspartate Aminotransferase2 assay are listed in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- Avoid hemolyzed samples due to high AST levels found in red blood cells. Specimens with hemoglobin levels greater than 10 mg/dL may cause falsely elevated results with the Aspartate Aminotransferase2 assay. Refer to the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert for further information.
- Specimens with triglyceride levels greater than 750 mg/dL may cause inaccurate results with the Aspartate Aminotransferase2 assay. Refer to the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert for further information.
- Potential interference has not been evaluated for substances other than those described in the SPECIFIC PERFORMANCE CHARACTERISTICS, Analytical Specificity, Interference section of this package insert.
- SmartWashes for assays impacted by Aspartate Aminotransferase2 must be configured to avoid interference due to reagent carryover. See the INSTRUMENT PROCEDURE section of this package insert for the required assay file updates.

### EXPECTED VALUES

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

#### Reference Range<sup>16</sup>

	Range (U/L)	Range (μkat/L)
Adults	11-34	0.18-0.57

### SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

Unless otherwise specified, the study results provided in this package insert were generated using the Calibration method.

## Precision

### Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A3.<sup>17</sup> Testing was conducted using 3 lots of the Aspartate Aminotransferase2 reagents, 3 lots of the Consolidated Chemistry Calibrator, 1 lot of commercially available controls, and 3 instruments. Two controls and 4 human serum panels were tested in a minimum of 2 replicates, twice per day, on 20 days on 3 reagent lot/calibrator lot/instrument combinations, where a unique reagent lot and a unique calibrator lot are paired with 1 instrument. The performance from a representative combination is shown in the following table.

Sample	n	Mean (U/L)	Within-Run (Repeatability)		Within-Laboratory <sup>a</sup>	
			SD	%CV	SD (Range <sup>b</sup> )	%CV (Range <sup>b</sup> )
Control Level 1	80	45	0.6	1.4	0.8 (0.8-1.3)	1.9 (1.9-3.0)
Control Level 2	80	228	0.9	0.4	1.6 (1.6-1.8)	0.7 (0.7-0.8)
Panel A	80	8	0.6	7.5	0.7 (0.7-1.0)	8.0 (8.0-12.1)
Panel B	80	33	0.7	2.1	0.9 (0.9-1.2)	2.8 (2.8-3.4)
Panel C	80	113	0.7	0.7	1.5 (1.5-2.2)	1.3 (1.3-1.9)
Panel D	80	3586	13.3	0.4	40.3 (40.3-47.3)	1.1 (1.1-1.3)

<sup>a</sup> Includes within-run, between-run, and between-day variability.

<sup>b</sup> Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot/instrument combinations.

Sample	n	Mean (μkat/L)	Within-Run (Repeatability)		Within-Laboratory <sup>a</sup>	
			SD	%CV	SD (Range <sup>b</sup> )	%CV (Range <sup>b</sup> )
Control Level 1	80	0.74	0.011	1.5	0.014 (0.014-0.023)	1.9 (1.9-3.0)
Control Level 2	80	3.80	0.014	0.4	0.026 (0.026-0.029)	0.7 (0.7-0.8)
Panel A	80	0.14	0.010	6.9	0.011 (0.011-0.017)	7.8 (7.8-12.1)
Panel B	80	0.54	0.012	2.2	0.015 (0.015-0.019)	2.8 (2.8-3.4)
Panel C	80	1.88	0.013	0.7	0.024 (0.024-0.036)	1.3 (1.3-1.9)
Panel D	80	59.78	0.221	0.4	0.672 (0.672-0.788)	1.1 (1.1-1.3)

<sup>a</sup> Includes within-run, between-run, and between-day variability.

<sup>b</sup> Minimum and maximum SD or %CV across the 3 reagent lot/calibrator lot/instrument combinations.

## Accuracy

A study was performed to estimate the bias of the Aspartate Aminotransferase2 assay relative to standard reference material (ERM-AD457/IFCC).

### Calibration Method

Testing was conducted using 3 lots of the Aspartate Aminotransferase2 reagents, 2 lots of the Consolidated Chemistry Calibrator, and 3 instruments. The bias ranged from -0.7 to 3.4% across all instruments, calibrator, and reagent lots.

### Calibration Factor Method

Testing was conducted using 3 lots of the Aspartate Aminotransferase2 reagents and 3 instruments. The bias ranged from -1.8% to 2.0% across all instruments and reagent lots.

## Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2.<sup>18</sup> Testing was conducted using 3 lots of the Aspartate Aminotransferase2 reagents on each of 2 instruments over a minimum of 3 days. The maximum observed limit of blank (LoB), limit of detection (LoD), and limit of quantitation (LoQ) values are summarized below.

	U/L	μkat/L
LoB <sup>a</sup>	2	0.03
LoD <sup>b</sup>	4	0.07
LoQ <sup>c</sup>	5	0.08

<sup>a</sup> The LoB represents the 95th percentile from  $n \geq 60$  replicates of zero-analyte samples.

<sup>b</sup> The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on  $n \geq 60$  replicates of low-analyte level samples.

<sup>c</sup> The LoQ is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met and was determined from  $n \geq 60$  replicates of low-analyte level samples.

## Linearity

A study was performed based on guidance from CLSI EP06-A.<sup>19</sup> This assay is linear across the analytical measuring interval of 5 to 4001 U/L (0.08 to 66.70 μkat/L).

## Analytical Specificity

### Interference

A study was performed based on guidance from CLSI EP07, 3rd ed.<sup>20</sup> Each substance was tested at 2 levels of the analyte (approximately 30 U/L and 120 U/L).

### Potentially Interfering Endogenous Substances

**No significant interference (interference within  $\pm 10\%$ )** was observed at the following concentrations.

No Significant Interference (Interference within $\pm 10\%$ )		
Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
Bilirubin - conjugated	40 mg/dL	474 μmol/L
Bilirubin - unconjugated	60 mg/dL	1026 μmol/L
Hemoglobin	10 mg/dL	0.10 g/L
Total protein	15 g/dL	150 g/L
Triglycerides	750 mg/dL	8.5 mmol/L

**Interference beyond  $\pm 10\%$  (based on 95% Confidence Interval [CI])** was observed at the concentrations shown below for the following substances.

Interference beyond $\pm 10\%$ (based on 95% Confidence Interval [CI])					
Potentially Interfering Substance	Interferent Level		Analyte Level		% Interference (95% CI)
	Default Units	Alternate Units	Default Units	Alternate Units	
Hemoglobin	15 mg/dL <sup>a</sup>	0.15 g/L	30 U/L	0.50 μkat/L	10% (8%, 12%)

<sup>a</sup> The hemoglobin interferent level presented in the table was generated using the Calibration Factor method.

### Potentially Interfering Exogenous Substances

**No significant interference (interference within  $\pm 10\%$ )** was observed at the following concentrations.

No Significant Interference (Interference within $\pm 10\%$ )		
Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
3-methyl-(triazene-1-yl)imidazole-4-carboxamide (MTIC)	0.6 mg/L	3.6 μmol/L
5-amino-4-imidazolecarboxamide (AIC)	3 mg/L	24 μmol/L
Acetaminophen	160 mg/L	1059 μmol/L
Acetylcysteine	150 mg/L	920 μmol/L
Acetylsalicylic acid	30 mg/L	167 μmol/L



No Significant Interference (Interference within $\pm 10\%$ )		
Potentially Interfering Substance	Interferent Level	
	Default Units	Alternate Units
Aminosalicilic acid (p-Aminosalicilic acid)	47 mg/dL	3074 $\mu\text{mol/L}$
Ampicillin-Na	80 mg/L	215 $\mu\text{mol/L}$
Ascorbic acid	60 mg/L	341 $\mu\text{mol/L}$
Biotin	4250 ng/mL	17 $\mu\text{mol/L}$
Ca-dobesilate	60 mg/L	143 $\mu\text{mol/L}$
Cefoxitin	6600 mg/L	15 mmol/L
Chlordiazepoxide	1 mg/dL	33 $\mu\text{mol/L}$
Clothiapine	9 mg/L	26 $\mu\text{mol/L}$
Cyclosporine	2 mg/L	1.7 $\mu\text{mol/L}$
Diclofenac	3 mg/dL	101 $\mu\text{mol/L}$
Doxycycline	20 mg/L	45 $\mu\text{mol/L}$
Furosemide	2 mg/dL	60 $\mu\text{mol/L}$
Hydroxocobalamin (Cyanokit)	500 mg/L	372 $\mu\text{mol/L}$
Ibuprofen	220 mg/L	1067 $\mu\text{mol/L}$
Levodopa	20 mg/dL	1014 $\mu\text{mol/L}$
Methylidopa	22 mg/dL	1041 $\mu\text{mol/L}$
Metronidazole	130 mg/L	759 $\mu\text{mol/L}$
Nitrofurantoin	0.3 mg/dL	13 $\mu\text{mol/L}$
Phenylbutazone	330 mg/L	1069 $\mu\text{mol/L}$
Pindolol	5 mg/L	20 $\mu\text{mol/L}$
Rifampicin	50 mg/L	61 $\mu\text{mol/L}$
Sodium Heparin	4 U/mL	N/A
Sulfapyridine	300 mg/L	1203 $\mu\text{mol/L}$
Sulfasalazine	300 mg/L	753 $\mu\text{mol/L}$
Suramin	500 mg/L	386 $\mu\text{mol/L}$
Temozolomide	5 mg/L	26 $\mu\text{mol/L}$
Theophylline (1,3-dimethylxanthine)	60 mg/L	333 $\mu\text{mol/L}$

N/A = Not applicable

**Interference beyond  $\pm 10\%$  (based on 95% Confidence Interval [CI])** was observed at the concentrations shown below for the following substances.

Interference beyond $\pm 10\%$ (based on 95% Confidence Interval [CI])					
Potentially Interfering Substance	Interferent Level		Analyte Level		%
	Default Units	Alternate Units	Default Units	Alternate Units	
Hydroxocobalamin (Cyanokit)	1000 mg/L	743 $\mu\text{mol/L}$	30 U/L	0.50 $\mu\text{kat/L}$	-12% (-14%, -11%)
Temozolomide	10 mg/L	52 $\mu\text{mol/L}$	30 U/L	0.50 $\mu\text{kat/L}$	12% (10%, 15%)

Interferences from medication or endogenous substances may affect results.<sup>21</sup>

## Method Comparison

A study was performed based on guidance from CLSI EP09-A3<sup>22</sup> using the Passing Bablok regression method.






Aspartate Aminotransferase2 on the ARCHITECT c System vs a comparator AST assay					
	n	Units	Correlation Coefficient	Intercept	Slope
Serum	168	U/L ( $\mu\text{kat/L}$ )	1.00	-1 (-0.01)	0.98
					Concentration Range 9-3910 (0.14-65.18)

## BIBLIOGRAPHY

- Burtis CA, Bruns DE, editors. *Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics*. 7th ed. St. Louis, MO: Elsevier;2015:322-323.
- Jacobs DS, De Mott WR, Oxley DK, editors. *Jacobs & DeMott Laboratory Test Handbook With Keyword Index*. 5th ed. Hudson (Cleveland), Ohio: Lexi-Comp; 2001:112-113.
- Kwo PY, Cohen SM, Lim JK. ACG clinical guideline: evaluation of abnormal liver chemistries. *Am J Gastroenterol* 2017;112(1):18-35.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Cuhadar S, Atay A, Koseoglu M, et al. Stability studies of common biochemical analytes in serum separator tubes with or without gel barrier subjected to various storage conditions. *Biochem Med* 2012;22(2):202-214.
- Heins M, Heil W, Withold W. Storage of serum or whole blood samples? Effects of time and temperature on 22 serum analytes. *Eur J Clin Chem Clin Biochem* 1995;33(4):231-238.
- Ikeda K, Ichihara K, Hashiguchi T, et al. Evaluation of the short-term stability of specimens for clinical laboratory testing. *Biopreservation & Biobanking* 2015;13(2):135-143.
- Cuhadar S, Koseoglu M, Atay A, et al. The effect of storage time and freeze-thaw cycles on the stability of serum samples. *Biochem Med* 2013;23(1):70-77.
- Clinical and Laboratory Standards Institute (CLSI). *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: CLSI; 2010.
- Westgard JO. *Basic QC Practices; Training in Statistical Quality Control for Medical Laboratories*. 4th ed. Westgard QC, Inc.; 2016.
- Schumann G, Bonora R, Ceriotti F, et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 5. Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. *Clin Chem Lab Med* 2002;40(7):725-733.
- Clinical and Laboratory Standards Institute (CLSI). *Establishing and Verifying an Extended Measuring Interval Through Specimen Dilution and Spiking; Approved Guideline—First Edition*. CLSI Guideline EP34. Wayne, PA: CLSI; 2018.
- Ceriotti F, Henny J, Queralto J, et al. Common reference intervals for aspartate aminotransferase (AST), alanine aminotransferase (ALT) and  $\gamma$ -glutamyl transferase (GGT) in serum: results from an IFCC multicenter study. *Clin Chem Lab Med*. 2010;48(11):1593-1601.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision of Quantitative Measurement Procedures: Approved Guideline—Third Edition*. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
- Clinical and Laboratory Standards Institute (CLSI). *Interference Testing in Clinical Chemistry*. 3rd ed. CLSI Guideline EP07. Wayne, PA: CLSI; 2018.
- Young DS. Laboratory test listings. In: *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. AACCPress; 2000:chap 3.
- Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

## ■ Key to Symbols

### ISO 15223 Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
<b>IVD</b>	<i>In Vitro</i> Diagnostic Medical Device
<b>LOT</b>	Lot Number
<b>REF</b>	List Number
<b>SN</b>	Serial number

### Other Symbols

<b>CONTAINS: AZIDE</b>	Contains Sodium Azide. Contact with acids liberates very toxic gas.
<b>DISTRIBUTED IN THE USA BY</b>	Distributed in the USA by
<b>FOR USE WITH</b>	Identifies products to be used together
<b>INFORMATION FOR USA ONLY</b>	Information needed for United States of America only
<b>PRODUCT OF IRELAND</b>	Product of Ireland
<b>R1</b>	Reagent 1
<b>R2</b>	Reagent 2
<b>Rx ONLY</b>	For use by or on the order of a physician only (applicable to USA classification only).

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

The ARCHITECT c System family of instruments consists of c4000, c8000, and c16000 instruments.

ARCHITECT and related brand marks are trademarks of Abbott.

Other trademarks are the property of their respective owners.



Abbott Ireland  
Diagnostics Division  
Lisnamuck, Longford  
Co. Longford  
Ireland  
+353-43-3331000



### **DISTRIBUTED IN THE USA BY**

Abbott Laboratories  
Abbott Park, IL 60064 USA

**Customer Service: Contact your local representative or find country-specific contact information on [www.corelaboratory.abbott](http://www.corelaboratory.abbott)**

For customers in the European Union: if, in the course of using this device, you have reason to believe that a serious incident has occurred, report it to the manufacturer and to your national authority.

Revised July 2021.

©2021 Abbott Laboratories