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
Port ID	Port Name	Order No.	Connector
1	LEAK TEST	H28611-1	
2	JET/AUX	H28611-26	
3	BIOPSY	H28611-3	
4	AIR	H28611-4	
5	ELEVATOR	H28611-5	
6	WATER	H28611-26	
7	SUCTION	H28611-7	
8	EXTRA	N/A	N/A
S	N/A	78400-930	

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
Hookup Connection Instructions for Use

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Hookup Diagram

Each Hookup is comprised of multiple ports and connector assemblies. Each port is labeled with a Port ID number (circled). This same number also identifies the corresponding Connector Assembly. Each Connector Assembly is fitted with a specific endoscope connector.

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Port ID	Port Name	Connector Assembly Description	Order No.
1	LEAK TEST	Leak Test Connector Assembly	H58112-1
2	JET/AUX	Air/Water Connector Assembly	N/A
3	BIOPSY	Biopsy Inlet Connector Assembly	H58112-3
4	AIR	Air Probe Connector Assembly	N/A
5	ELEVATOR	Auxiliary Water Connector Assembly	N/A
6	WATER	Air/Water Connector Assembly	N/A
7	SUCTION	Biopsy 2 Inlet Connector Assembly	H58112-7
8	EXTRA	N/A	N/A
S	N/A	N/A	N/A

Connector Assemblies are specific to individual hookups.

Designated Hookup

For the Following PENTAX® Endoscopes*					
CHF-10	CHF-B260	CHF-CB20	CHF-CB30L	CHF-CB30S	CHF-P10
CYF	CYF-2	CYF-200	CYF-240	CYF-3	CYF-4
CYF-5	CYF-5R	CYF-V	CYF-V2	CYF-V2R	CYF-VH
CYF-VHR	HYF-P	HYF-V	HYF-XP	SIF-SW	URF-P
URF-P2	URF-P3	URF-P5	URF-V	URF-P7	URF-P7R
URF-P2	URF-P3	URF-P5	URF-V	URF-P7	URF-P7R
URF-V3	URF-V3R				

*This list of corresponding endoscope models is subject to change without notice. For a full list, please go to: www.medivators.com/hookuplookup

*Not all endoscopes may be available in every region. Contact your local representative for device availability.

Intended Purpose

Hookups are intended to provide the interface between an Automatic Endoscope Reprocessor (AER) and endoscopes, which allow perfusion of liquids dispensed by the AER to contact and high-level disinfect the channels of flexible endoscopes.

Clinical Benefit

This ability of the device to flush endoscope channels with high-level disinfection during the reprocessing cycle of fully immersible endoscopes helps reduce post-endoscopic infections, thereby ensuring patient safety.

Warnings & Precautions

- Ensure proper alignment of all connectors with the endoscope mating parts prior to attachment.
- This hookup is intended for connection only with the specific endoscope models identified in the appropriate Medivators online hookup guide. Use of this hookup with endoscopes other than those specified by the hookup application guide may result in inadequate HLD.
- Users operating Medivators AERs and hookups must be trained and competent in the understanding of endoscope channel systems.
- Prior to reprocessing, users must verify that the hookup used contains connections for all appropriate channels that require a separate reprocessing connection according to the hookup instructions.
- All connectors and adapters must remain firmly attached and unrestricted for the entire disinfection cycle to ensure adequate HLD. Users must inspect the endoscope and hookup combination to verify proper connection and flow.
- Modifications or repairs to hookups which do not correspond to the manufacturer's specifications may result in inadequate HLD and/or damage to endoscopes.
- Failure to properly clean and prepare endoscope for immersion prior to hookup connection may result in inadequate HLD and/or damage to the endoscope.
- Endoscopes must be inspected for damage and verified to be in proper working order prior to hookup connection. Connection of hookups to damaged endoscopes or endoscopes with obstructed/restricted channels could result in inadequate HLD.
- Manual leak testing must be carried out in accordance with manufacturer's recommendations prior to reprocessing.
- Install all necessary waterproof caps, plugs and cleaning adapters prior to immersion.
- Do not autoclave Medivators hookups.
- Medivators makes no claim on high-level disinfection efficacy when these reprocessing instructions are not followed, or when this hookup kit is applied to endoscopes other than those specified in the online interactive hookup guide.
- Avoid biological contamination. Always wear personal protective equipment when handling endoscopes.
- An Advantage Hookup Cassette combines an endoscope tray with the hookup connector block and are designated with a CAS suffix (applicable to metal cassette only). Specific endoscope connections for this hookup are identical to non-cassette versions.

Adverse Event Reporting

User facilities shall report adverse events to Medivators, a Cantel Medical company, as well as the competent authority of the EU member state in which the user is established, if applicable.

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
For additional symbols reference information, definition and translation, please reference the Cantel Symbols Glossary on the Medivators website, www.medivators.com.

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REF 5-8-112, 5-8-112HAN, 5-8-112CAS

Locate the endoscope port and follow the instructions for attaching the corresponding Endoscope Connector. All connectors must be attached. Ensure the hookup block is appropriately connected to the base of the AER.

Biopsy Connector
Align stainless steel fitting with biopsy port, push and turn clockwise to connect.

Leak Test Connector
Align the slot on connector with the pin on the scope. Push on and turn clockwise to engage. To remove leak test connector, turn connector, counter-clockwise until slot is aligned with pin and pull off. Do not pull until slot is aligned with pin.

This hookup is designed to be used with two identical endoscope models. Channel connections are duplicates of each other. This hookup cannot be used with only one endoscope or dissimilar models.

Representative endoscope for illustrative purposes only. For exact port locations, refer to endoscope manufacturer's instruction guide.

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Catalog number

Batch code

Authorized representative in the European Community

Consult Instructions for Use or consult electronic Instructions for Use

Manufacturer

Date of manufacture

Importer

Medical device

Unique device identifier

ID

Australian sponsor

Lock and unlock

Leak test

Jet/Aux Fw Jet

Biopsy

Air

Elevator Fw Jet

Water

Suction

Extra


See AER symbols glossary or user manual for more information.

Hookup Connector Assembly Replacement Instructions

Using MEDIVATORS™ ADVANTAGE PLUS™ Reprocessors, ADVANTAGE PLUS™ Pass-Thru Reprocessors and ISA™ Reprocessors.

REF 2-8-332 / 2-8-332HAN / 2-8-332CAS

Hookup Diagram

Directions for Replacement of Connector Assemblies

1. Identify if the Hookup Connector Block is the new (Figure 1 below) or previous (Figure 2 below) design.
2. Identify current Connector Assembly on the Hookup Connector Block to be replaced.
3. Remove current Connector Assembly by pulling associated tubing **B** from the Hookup Connector Block Barb **A** (NOTE: If using the previous design, unscrew the Connector Coupling **F** of Connector Assembly from the Port Connector **C**, keeping the Port Connector firmly held in place and pull the tubing off the Port Connector stem **D**). Discard this Connector Assembly.
4. Retrieve the enclosed new Connector Assembly. NOTE: that it may not be identical to the Connector Assembly it is replacing. Discard Ferrule **E** and Connector Coupling **F** if using the new design. Confirm Connector Assembly labeling matches the Hookup Connector Block Label Port Name and the Port ID Number using the Connector Assembly Identification Chart. NOTE: Port ID Number may not be identified on some older Hookup Connector Block labels.
5. Push tubing **B** of the new Connector Assembly onto the Hookup Connector Block Barb **A**. NOTE: If experiencing difficulty in attaching the tubing to the Hookup Connector Block Barb **A**, dip the tip of the tubing into hot water for 10 seconds. (NOTE: If using the previous design, push tubing of new Connector Assembly onto the Port Connector Stem **C** and manually screw Connector Coupling **F** onto Port Connector **C** until tight. Do not over tighten. Ensure Port Connector **C** does not turn. If needed, use a wrench to prevent turning.)

Connector Assembly Identification Chart

Port ID	Port Name	Order No.	Connector
1	LEAK TEST	H28332-1	
2	JET/AUX	N/A	N/A
3	BIOPSY	H28332-3	
4	AIR	N/A	N/A
5	ELEVATOR	N/A	N/A
6	WATER	H28332-6	
7	SUCTION	H28332-7	
8	EXTRA	N/A	N/A
S	N/A	N/A	N/A

If a Connector Assembly needs replacement, locate the Port Number/Name and use the corresponding order number found in the chart above. Connector Assemblies are specific to individual hookups.

Hookup Connector Block Types

Exploded view of Port Connector, Coupling Connector, Ferrule and Tubing:

Figure 1 - New Design



Figure 2 - Previous Design

NOTE: Ferrule must be present to prevent leaking.

Warnings & Precautions

- This hookup is intended for connection only with the specific endoscope models identified in the appropriate Medivators online hookup guide. Use of this hookup with endoscopes other than those specified by the online hookup guide may result in inadequate HLD.
- Users operating Medivators AERs and hookups must be trained and competent in the understanding of endoscope channel systems.
- Prior to reprocessing, users must verify that the hookup used contains connections for all appropriate channels that require a separate reprocessing connection according to the hookup instructions.
- All connectors and adapters must remain firmly attached and unrestricted for the entire disinfection cycle to ensure adequate HLD. Users must inspect the endoscope and hookup combination to verify proper connection and flow.
- Modifications or repairs to hookups which do not correspond to the manufacturer's specifications may result in inadequate HLD and/or damage to endoscopes.
- Failure to properly clean and prepare endoscope for immersion prior to hookup connection may result in inadequate HLD and/or damage to the endoscope.
- Endoscopes must be inspected for damage and verified to be in proper working order prior to hookup connection. Connection of hookups to damaged endoscopes or endoscopes with obstructed/restricted channels could result in inadequate HLD.
- Manual leak testing must be carried out in accordance with manufacturer's recommendations prior to reprocessing.
- Install all necessary waterproof caps, plugs and cleaning adapters prior to immersion.
- Do NOT autoclave Medivators hookups.
- Medivators makes no claim on high-level disinfection efficacy when these reprocessing instructions are not followed, or when this hookup kit is applied to endoscopes other than those specified in the online interactive hookup guide.
- Avoid biological contamination. Always wear personal protective equipment when handling endoscopes.

If you have any questions regarding this hookup, in the United States, please contact Cantel Medivators Technical Support Representative at 1-800-444-4729. Outside the United States, please contact your Medivators representative.

ADVANTAGE PLUS™, ISA™ and MEDIVATORS™ are trademarks of Cantel Medical Corp., its affiliates or related companies. All other company or product names referenced are trademarks of their respective owners.

 Medivators Inc.
14605 28th Avenue North
Minneapolis, MN
55447-4822 USA
Tel: +1.763.553.3300
Toll Free: +1.800.444.4729
Fax: +1.763.553.3387
www.medivators.com

 Cantel Medical (Italy) S.r.l.
Via Laurentina, 169
00071 Pomezia (RM) Italy
Tel: +39.06.9145399
Fax: +39.06.9146099
www.cantelmedical.eu

Cantel (Australia) Pty Ltd 
9 Arco Lane, Heatherton,
Victoria 3202, Australia
Toll Free: +61.1300.211.422
Fax: +61.3.9088.8342
Email: info@cantel.com.au 

THE INTERNATIONAL CERTIFICATION NETWORK

CERTIFICATE

CISQ/IMQ has issued an IQNet recognized certificate that the organization:

CANTEL MEDICAL (ITALY) SRL

VIA LAURENTINA 169 - 00071 POMEZIA (RM)

has implemented and maintains a

Quality Management System

for the following scope:

Design, development, manufacturing of disinfectants, sterilizers and detergents for medical devices. Design, development, production, sales and technical service of device for washing, disinfection and sterilization of medical devices. Design, development, production management of conservation and transport systems for endoscopes

Further clarifications regarding the applicability of UNI CEI EN ISO 13485:2016 requirements may be obtained by consulting the organization

which fulfills the requirements of the following standard:

UNI CEI EN ISO 13485:2016

Issued on: 2021 - 01 - 21

Expires on: 2024 - 07 - 05

This attestation is directly linked to the IQNet Partner's original certificate and shall not be used as a stand-alone document

Registration Number: IT - 126041

Alex Stoichitoiu
President of IQNET

Ing. Mario Romersi
President of CISQ

IQNet Partners*:

AENOR Spain AFNOR Certification France APCER Portugal CCC Cyprus CISQ Italy
CQC China CQM China CQS Czech Republic Cro Cert Croatia DQS Holding GmbH Germany EAGLE Certification Group USA
FCAV Brazil FONDONORMA Venezuela ICONTEC Colombia Inspecta Sertifiointi Oy Finland INTECO Costa Rica
IRAM Argentina JQA Japan KFQ Korea MIRTEC Greece MSZT Hungary Nemko AS Norway NSAI Ireland
NYCE-SIGE México PCBC Poland Quality Austria Austria RR Russia SII Israel SIQ Slovenia
SIRIM QAS International Malaysia SQS Switzerland SRAC Romania TEST St Petersburg Russia TSE Turkey YUQS Serbia

www.imq.it

IQNet, the association of the world's first class certification bodies, is the largest provider of management System Certification in the world. IQNet is composed of more than 30 bodies and counts over 150 subsidiaries all over the globe.

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UNITA' OPERATIVE / OPERATIVE UNITS

VIA LAURENTINA 169 - 00071 POMEZIA (RM)

E' CONFORME ALLA NORMA / IS IN COMPLIANCE WITH THE STANDARD

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Certificate No. 14038-9-2020

CERTIFICATE TO FOREIGN GOVERNMENT

In order to allow the importation of United States products into foreign countries, the U.S. Food and Drug Administration (FDA) certifies the following information concerning the product(s) to be exported listed below:

Name of Product(s)

See Attached List

(Three Pages)

Name of Manufacturer/Distributor, Address

See Attached List

(One Page)

The product(s) described above (and the manufacturing/distribution site(s) which produces/distributes it) is subject to the jurisdiction of the FDA under the Federal Food, Drug, and Cosmetic Act.

It is certified that the above product(s) may be marketed in, and legally exported from, the United States of America at this time. The manufacturing plant(s) in which the product(s) is produced is subject to periodic inspections. The last such inspection showed that the plant(s), at that time, appeared to be in substantial compliance with current good manufacturing practice requirements for the product(s) listed above.

Sincerely,

CDR Cesar A. Perez, PhD, Director
DRP2: Division of Establishment Support
Office of Regulatory Programs
Office of Product Evaluation and Quality
Center for Devices and Radiological Health
U.S. Food and Drug Administration, DHHS

This certificate is valid from September 18, 2020 to September 17, 2022.

Certificate No. 14038-9-2020

Certificate to Foreign Government - Name of Manufacturer/Distributor Attachment Page 1 of 1

Name of Owner Operator

MEDIVATORS
14605 28TH AVE NORTH
Minneapolis, MN USA 55447

Name of Manufacturer

Medivators Inc
14605 28th Ave N
MINNEAPOLIS, MN USA 55447

Medivators Inc
3150 Pollok Drive
Conroe, TX USA 77303

----END OF MANUFACTURER/DISTRIBUTOR LIST----

Certificate No. 14038-9-2020

Certificate to Foreign Government - Name of Product(s) Attachment Page 1 of 3

Name of Owner Operator

MEDIVATORS
14605 28TH AVE NORTH
Minneapolis, MN USA 55447

Name of Manufacturer

Medivators Inc
14605 28th Ave N
MINNEAPOLIS, MN USA 55447

Medivators Inc
3150 Pollok Drive
Conroe, TX USA 77303

Name of Product(s)

Renatron PA Dialyzer Reprocessing System
Renatron PA 100 Dialyzer Reprocessing System-- RS-8315
Renatron II Dialyzer Reprocessing System
Renatron II 100 Dialyzer Reprocessing System-- RS-8335
Renalog RM-- 78398-192, 78398-204
Renalin Cold Sterilant-- 78397-521
Renalin PA Cold Sterilant-- 78397-638, 78397-639, 78397-452
Renalin 100 Cold Sterilant-- 78398-060, 78397-844, 78397-845, 78397-846, 78397-970
Renalin 100 Cold Sterilant with test kit-- 78397-845, 78397-846, 78397-970
Renalin Indicator Test Strips --78199-000
Renalin Residual Test Strips-- 78198-000
Perassay 500 Peracetic Acid Test Strips-- 78378-000
Actril Cold Sterilant-- 78337-000, 78399-667
Actril Cold Sterilant with test kit-- 78270-000
Actril Residual Test Strips-- 78258-000
Actril Indicator Test Strips-- 78259-000
Hemocor HPH Hemoconcentrator; HPH Mini, HPH Junior
Hemocor HPH Hemoconcentrator; HPH 400, HPH 400TS
Hemocor HPH Hemoconcentrator; HPH 700, HPH 700TS
Hemocor HPH Hemoconcentrator; HPH 1000, HPH 1000TS
Hemocor HPH Hemoconcentrator; HPH 1400, HPH 1400TS
Renafo II Hemofilter; HF Minifilter Plus
Renafo II Hemofilter; HF Junior
Renafo II Hemofilter; HF400
Renafo II Hemofilter; HF700
Renafo II Hemofilter; HF1200
Renafo II Hemofilter; HF2000
Primus Hollow Fiber Dialyzer; Primus 1350
Fiberflo Hollow Fiber Cartridge Water Filter--FF50, FF100, FF200
Bicarbonate Concentrate Liquid; BC-1-L Renasol, MB-330-L Centrisol
Bicarbonate Concentrate Powder; Renasol BC-1, BC-1-15, BC-1-25
Bicarbonate Concentrate Powder; Centrisol MB-330, MB-330-15, MB-330-25
Acetate Concentrate Nephrosol; Acetate Concentrate
Renasol Acid Concentrate - SB-1003, SB-1014, SB-1019, SB-1020, SB-1030, SB-1057, SB-1063, SB-1078, SB-1086
Centrisol Acid Concentrate - SB-111, SB-119, SB-123, SB-127, SB-129, SB-140, SB-142, SB-143, SB-144, SB-145
Medivators Advantage Plus Endoscope Reprocessing System
Medivators Advantage Automated Endoscope Reprocessor
MDS Modular Disinfection System
ADVANTAGE PLUS Pass-Thru Automatic Endoscope Reprocessor -- ADVPT-3005, ADVPT-3006
DSD Series Automatic Endoscope Reprocessor-- DSD-201(LT)(Model No.: DSD-2007, DSD-2011), SSD-102(LT)

Certificate No. 14038-9-2020

Certificate to Foreign Government - Name of Product(s) Attachment Page 2 of 3

DSD Series Automatic Endoscope Reprocessor-- DSD Edge
CER Series Endoscope Reprocessor-- CER-1, CER-2, CER-1 Optima, CER-2 Optima
Veriscan Automated Endoscope Leak Detection System-- VER-LT
Rapicide PA High-Level Disinfectant-- ML02-0117, ML02-0116, ML02-0115
Rapicide OPA-28 High-Level Disinfectant-- ML02-0127, ML02-0132, ML02-0138
Rapicide PA RTU High-Level Disinfectant & Sterilant- ML02-0125
Rapicide High-Level Disinfectant and Sterilant-- ML02-0059, ML02-0108
Rapicide PA Disinfectant & Sterilant Test Strips - ML02-0118
Rapicide PA RTU High-Level Disinfectant & Sterilant Test Strips - ML02-0129
Rapicide OPA-28 Test Strips-- ML02-0137, ML02-0136
Rapicide Test Strips-- ML02-0120
MediClean EZ-- ML02-0131
Intercept Detergent-- ML02-0106, ML02-0134
Intercept Plus - ML02-0145
Intercept Wipes-- ML02-0107
Intercept Foam - ML02-0139
Scope Buddy Endoscope Flushing Aid-- EFA-US-G, EFA-IN-G (Model No. ECA-100G)
Scope Buddy Tubing Accessory Kit- 78399-669
Minnicare Cold Sterilant -78398-229, 78397-983, 78397-825
Minnicare Residual Test Strips -78338-000
Minnicare Peracetic Acid Test Strips- 78339-000
EndoStratus Irrigation Pump EGA-500, EGA-500E, EGA-500P
EndoGator Irrigation Pump EGP-100, EGP-100E, EGP-100P
EndoStratus CO2 Insufflator EGA-501, EGA-501E, EGA-501P
EndoStratus Irrigation Pump and CO2 Insufflator System -- EGA-SYST
EndoSmartCap Irrigation Tubing and Accessories 100145, 100145U, 100145CO2, 100145CO2U
ENDO SMARTCAP Co2 Tubing 100145CO2EX, 100150CO2EX, 100165CO2EX
EndoSmartCap Irrigation Tubing and Accessories 100150CO2P
EndoSmartCap Irrigation Tubing and Accessories 100150, 100150U, 100150CO2, 100150CO2U
EndoSmartCap Irrigation Tubing and Accessories 100160, 100160U, 100160P, 100160S
EndoSmartCap Irrigation Tubing and Accessories 100165, 100165U, 100165CO2, 100165CO2U
EndoSmartCap Irrigation Tubing and Accessories 100140, 100141, 100551, 100551P
EndoSmartCap Irrigation Tubing and Accessories 100162
EndoSmartCap Adapter 100555
ENDOGATOR Cartridge 100110
EndoGator Hybrid Irrigation Tubing 100605, 100605U, 100605S, 100606, 100606U, 100606S
EndoGator Hybrid Irrigation Tubing 100608, 100608U, 100608S, 100609, 100609U, 100609S
EndoGator Hybrid Irrigation Tubing 100610, 100610U, 100610S
EndoGator Hybrid Irrigation Tubing 100630, 1001630U, 100631, 100631U
EndoGator Irrigation Tubing and Accessories Kits 100130P, 100130, 100130U, 100130F, 100130FU
EndoGator Irrigation Tubing and Accessories Kits 100125, 100125F, 100126
EndoGator Irrigation Tubing and Accessories Kits 200230P
EndoGator Irrigation Tubing and Accessories Kits 200230, 200230U, 200230F, 200230FU
EndoGator Irrigation Tubing and Accessories Kits 100135, 100136, 100105
EndoGator Irrigation Tubing and Accessories Kits 100131, 100115, 100116, 100116P
EndoGator Irrigation Tubing and Accessories Kits 100241, 100242, 100242P, 100115, 100116
EndoGator Irrigation Tubing and Accessories Kits 100600, 100601, 100602
EndoGator Irrigation Tubing and Accessories Kits 100603, 100611, 100612, 100613, 100614
EndoGator Irrigation Tubing and Accessories Kits 100615, 100616, 100618
EndoGator Irrigation Tubing and Accessories Kits 100620, 100621, 100622, 100623
EndoGator Irrigation Tubing and Accessories Kits 100625, 100626, 100627, 100628, 100629
EndoGator Irrigation Tubing and Accessories Kits 100630, 100630S, 100631, 100631S
Defendo Disposable Biopsy, Air/Water and Suction Valves 100301, 100302, 100303, 100305
Defendo Disposable Biopsy, Air/Water and Suction Valves 100306, 100310, 100312, 100313
Defendo Disposable Biopsy, Air/Water and Suction Valves 100314, 100315, 100316, 100317
DEFENDO Valve and Connector Kit for Olympus GI Endoscopes 100311
DEFENDO biopsy valve 100318, 100319, 100320
Defendo Olympus Kit 100322, 100323

Certificate No. 14038-9-2020

Certificate to Foreign Government - Name of Product(s) Attachment Page 3 of 3

Endo Carry-on Kit 100800, 100801, 100802, 100803, 100305, 100310

Endoscope Cleaning Brushes and Pull Thru Channel Cleaning Devices 100401, 100402, 100403

Endoscope Cleaning Brushes and Pull Thru Channel Cleaning Devices 100405, 100406

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
A handwritten signature in black ink, appearing to be 'S. O.', written over a horizontal line.

DocuSign

IMQ

EC CERTIFICATE

Certificate No 1812/MDD

Full Quality Assurance System Approval Certificate

On the basis of our examination carried out according to Annex II, excluding section 4, of the Directive 93/42/EEC and its revised version, we hereby certify that:

CANTEL MEDICAL (ITALY) SRL

00071 POMEZIA (RM) - VIA LAURENTINA 169 (ITA) - Italy

manages in the factory of:

00071 POMEZIA (RM) - VIA LAURENTINA 169 (ITA) - Italy

a quality assurance system ensuring the conformity of the following products:

Cold chemical washer disinfectant and sterilizer for endoscopes

Cold chemical sterilant for medical devices

Disinfectants for medical devices

Multi-enzyme detergent, decontaminant disinfectant for medical devices

Disinfectants, decontaminants and detergents for medical devices

Disinfectants and detergents for medical devices

Decontaminants and disinfectants for medical devices

Storage and transport systems for endoscopes

Washer disinfectant for endoscopes

series and type refs in the Annex

with the relevant essential requirements of the aforementioned directive (from design to final inspection and testing) and it is subject to surveillance as specified in section 5 of Annex II. For class III devices, this certificate is valid only with the relevant EC Design-Examination Certificate of Annex II.4.

Reference to IMQ files Nos:

DM15A0449933-01; DM15E0572628-01; DM16A0607476-01; DM16-0000589; DM16-0002190-01; DM19-0043082-01; DM19-0043104-01; DM20-0050214-01; DM20-0048482-01; DM20-0051086-01; DM20-0047938-01.

This Approval Certificate is issued by IMQ S.p.A. as Notified Body for the Directive 93/42/EEC and its revised version. Notified Body notified to European Commission under number: 0051.

Date: 2015-07-20
 Updated: 2020-05-08
 Substitution Date: 2020-04-07
 Expiry Date: 2024-05-26

DocuSign

IMQ

EC CERTIFICATE

Certificate No 1812/MDD

Annex

Cold chemical washer disinfectant and sterilizer for endoscopes

Type ref. MEDIVATORS ISA
Trade mark Cantel Medical (Italy) S.r.l.

Cold chemical sterilant for medical devices

Type ref. ADASPOR PRONTO; ADASPOR CONCENTRATO, ADASPOR M CONCENTRATO; ADASPOR MONODIE; ADASPOR PENTADIE; ADASPOR SINGLE SHOT; PROLYSTICA AUTO PAA; ISASPOR SINGLE SHOT; ADASPOR PLUS SINGLE SHOT, ADASPOR PLUS PRONTO (READY TO USE); ADASPOR PLUS CONCENTRATO; ADASPOR PLUS MONODIE; ADASPOR PLUS PENTADIE, ADASPOR PLUS M CONCENTRATO.
Trade mark CANTEL

Disinfectants for medical devices

Type ref. BLUESTERIL ALCOLICO; BLUESTERIL FERRI; BLUESTERIL SPRAY.
Trade mark CANTEL

Multi-enzyme detergent, decontaminant disinfectant for medical devices

Type ref. NEO PROTEOZIM PLUS 500; PROTEOZIM PLUS 400.
Trade mark CANTEL

Disinfectants, decontaminants and detergents for medical devices

Type ref. ISACLEAN, PROTEODONT.
Trade mark CANTEL

Disinfectants and detergents for medical devices

Type ref. BACTRYL SPRAY; BACTRYL WIPES; ISACLEAN SPRAY; SPOREXIN SPRAY; SPOREXIN WIPES; SPOREXIN VACUUM.
Trade mark CANTEL

Decontaminants and disinfectants for medical devices

Type ref. PROTEAZONE; PROTEAZONE OD.
Trade mark CANTEL

Storage and transport systems for endoscopes

Type ref. CLEANASCOPE; CLEANASCOPE ADVANTAGE.
Trade mark CANTEL

·
·
·

Date: 2015-07-20
Updated: 2020-05-08
Substitution Date: 2020-04-07
Expiry Date: 2024-05-26

DocuSign
IMQ

EC CERTIFICATE

Certificate No 1812/MDD

Annex

Washer disinfector for endoscopes

Type ref. INNOVA E3s; INNOVA E3s CMS; INNOVA E4s CMS.
Trade mark CANTEL

Date: 2015-07-20
Updated: 2020-05-08
Substitution Date: 2020-04-07
Expiry Date: 2024-05-26

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The undersigned CANTEL MEDICAL (ITALY) S.R.L.. declares that the medical device **MEDIVATORS ISA** meets the requirements of Annex II of EEC Directive 93/42 and its revised version and meets the applicable provisions thereof with the relevant essential requirements of the aforementioned directive from design to final inspection and testing.

CANTEL MEDICAL (ITALY) S.R.L. also claims that the product meets the following standards:


UNI EN ISO 15883-1
UNI EN ISO 15883-4
UNI CEN ISO/TS 15883-5
IEC 61010-1
IEC 61010-2-040
EN 61326-1
CEI EN 62366


Data/date 11/09/2017

GIORNO/DAY - MESE/MONTH - ANNO/YEAR

MANAGING DIRECTOR

V. M. SUPEKAR

NSAI

Certificate of Registration of Quality Management System to I.S. EN ISO 13485:2016

The National Standards Authority of Ireland certifies that:

Medivators Inc.
3150 Pollok Dr.
Conroe, TX 77303
USA

has been assessed and deemed to comply with the requirements of the above standard in respect of the scope of operations given below:

The Design and Development, Manufacture, and Distribution of Medical Devices used in Endoscopy Procedures and General Surgery including Disposable Valves, Sterile Tubing Sets, Endoscopy Cuffs, Endoscopy Procedure Kits, and related accessories.

Additional sites covered under this multi-site certification are listed on the Annex (File No. MD19.4716)

Approved by:
Geraldine Larkin
Chief Executive Officer

Approved by:
Caroline Dore Geraghty
Director of Medical Devices
/ Head of Notified Body

Registration Number: MD19.4716
Certification Granted: May 30, 2012
Effective Date: June 02, 2020
Expiry Date: May 15, 2022

NSAI

Annex to Certificate Number: MD19.4716

Scope of Registration:

The Design and Development, Manufacture, and Distribution of Medical Devices used in Endoscopy Procedures and General Surgery including Disposable Valves, Sterile Tubing Sets, Endoscopy Cuffs, Endoscopy Procedure Kits, and related accessories.

Activity

Location

HQ, Administration, Design & Development, Manufacturing, Distribution

Medivators Inc.
3150 Pollok Dr.
Conroe, TX 77303
USA
File No.: MD19.4716

Design & Development


Medivators Inc.
3101 Pollok Dr.
Conroe, TX 77303
USA
File No.: MD19.4716/A

Administration, Manufacturing

Medivators Inc.
501 North Farm to Market 3083 Rd E
Conroe, TX 77303
USA
File No.: MD19.4716/B

Manufacturing, Distribution

Medivators Inc.
750 Conroe Park North
Conroe, TX 77303
USA
File No: MD19.4716/C

NSAI

Annex to Certificate Number: MD19.4716

Scope of Registration:

The Design and Development, Manufacture, and Distribution of Medical Devices used in Endoscopy Procedures and General Surgery including Disposable Valves, Sterile Tubing Sets, Endoscopy Cuffs, Endoscopy Procedure Kits, and related accessories.

Activity

Location

Distribution

Medivators Inc.
3155 Conroe Park North
Conroe, TX 77303
USA
File No: MD19.4716/D

**Verified by:
Operations Manager**