LISS/Coombs

ID-Card

B004014 02.13

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Indirect and direct antiglobulin test

Product-Identification: 50531

INTRODUCTION

Polyspecific anti-human globulin (AHG) reagents are used for routine alloantibody detection and identification, compatibility tests and the direct antiglobulin test (DAT).

The most important function of the polyspecific AHG reagent is to detect the presence of IgG. The importance of anticomplement in the AHG reagent is debatable since antibodies detectable only by their ability to bind complement are rather rare. However, anti-C3d activity is important for the DAT in the investigation of autoimmune hemolytic anemia (AIHA) [4]. A positive DAT generally indicates that the red cells are coated *in vivo* with immunoglobulin and/or complement.

The microtubes of the ID-Card "LISS/Coombs" contain polyspecific AHG, to be used for antibody screening, antibody identification, crossmatch and the DAT. For the indirect antiglobulin test (IAT), labour intensive washing procedures are eliminated, due to the fact that the red cell suspension is added to the microtube before the plasma/serum, creating a barrier over the gel suspension, thus avoiding neutralization of the AHG by serum IgG proteins.

The anti-human globulin IgG used in the ID-Card "LISS/Coombs" is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

The ID-Card "LISS/Coombs" is suitable for the DAT, for the compatibility test, for antibody screening and identification with "ID-DiaCell" and "ID-DiaPanel".

REAGENTS**IVD**

ID-Card "LISS/Coombs" with 6 microtubes containing polyspecific AHG (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9) within the gel matrix.

Preservative: < 0.1% NaN₃.

Do not store near any heat, air conditioning sources or ventilation outlets.

Stability: see expiry date on label.

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspension.
- ID-DiaCell, ID-DiaPanel: Test cell reagents.

(see related package insert)

FURTHER MATERIALS REQUIRED

- ID-Dispenser
- ID-Pipetor
- ID-Tips (pipetor tips)
- Suspension Tubes
- ID-Working table
- ID-Incubator 37 °C
- ID-Centrifuge 6, 12 or 24

SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.

When the use of serum instead of plasma is required, the serum must be well cleared, by centrifugation at 1500 g for 10 minutes, before use avoid fibrin residues, which may interfere with the reaction pattern.

PREPARATION OF BLOOD SAMPLE

a) Red cell suspension for DAT or autocontrol

Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:

Allow the diluent to reach room temperature before use.

1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl of packed red cells, mix gently.

The cell suspension may be used immediately.

b) Red cell suspension for crossmatch procedures

Prepare a 0.8% red cell suspension in ID-Diluent 2 as above. Where whole blood directly from a segment of the blood bag is used, add 20 µl of blood to the 1.0 ml of ID-Diluent 2.

c) Plasma or serum for indirect antiglobulin test (IAT) procedures

Where samples are not for immediate testing they should be stored at 2–8 °C after separation (see also under "Sample material") for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

TEST PROCEDURES

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.

1. Direct antiglobulin test (DAT)

1. Identify the appropriate microtubes of the ID-Card "LISS/Coombs" with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.

3. Pipette 50 µl of the red cell suspension to the appropriate microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

II. Antibody screening (IAT)

Use ready-to-use test cell reagents "ID-DiaCell".

Allow the test cell reagents and samples to reach room temperature before use.

1. Identify the appropriate microtubes of the ID-Card "LISS/Coombs" with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each test cell reagent to the appropriate microtubes (marked with the corresponding test cell).
4. When an autocontrol is to be included, pipette 50 µl of the sample's own red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's or donor's plasma or serum to each microtube.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

III. Antibody Identification (IAT)

Use the ready-to-use test cell reagent "ID-DiaPanel".

Allow the test cell reagents and samples to reach room temperature before use.

1. Identify two ID-Cards "LISS/Coombs" with the patient's or donor's name or number.
 2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
 3. Pipette 50 µl of each "ID-DiaPanel" test cell to the appropriate microtube (marked 1 to 11).
 4. Pipette 50 µl of the sample's own red cell suspension to the 12th microtube (autocontrol).
-

5. Add 25 µl of the patient's or donor's plasma or serum to all 12 microtubes.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

IV. Compatibility test

1. Identify the appropriate microtubes of the ID-Card "LISS/Coombs" with recipient's and donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the donor red cell suspensions to the appropriate microtubes.
4. For the autocontrol, pipette 50 µl of the patient's own red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's plasma or serum to each microtube.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

INTERPRETATION OF THE RESULTS

A) Principle

Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.

Negative: Compact button of cells on the bottom of the microtube.

B) Reactions for:

1. Direct antiglobulin test (DAT)

- A negative reaction indicates absence of detectable IgG antibodies or C3d complement component on the red cells.
- A positive reaction (± to +++) indicates that the patient's red cells are sensitized (red cells coated with IgG antibodies and/or C3d).

II. Antibody screening

- A negative reaction indicates the absence of detectable irregular antibodies in the patient's or donor's serum or plasma.
- A positive reaction indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents "ID-DiaCell I-II" or "ID-DiaCell I-II-III" corresponds to the lot number indicated on the antigen table.
- Following the reaction pattern and the antigen configuration, the type of antibody present may be indicated. Perform the usual further tests to identify the antibody.
- A positive reaction with one or more test cells and a negative autocontrol suggest the presence of a specific antibody (see Remarks).
- A positive reaction with all test cells and a positive autocontrol may be due to non-specific reactions.
- A positive reaction with all test cells and a positive autocontrol but with one or more test cells showing a stronger positive reaction than the autocontrol, the patient sample should be submitted for further testing, to investigate the possibility of an underlying allo-antibody.

III. Antibody Identification

- A positive reaction indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents "ID-DiaPanel" corresponds to the lot number indicated on the antigen table.
 - Following the reaction pattern and the antigen configuration, the type of antibody present can, in most cases, be identified (autocontrol must be negative).
 - A positive reaction with all "ID-DiaPanel" test cells and a negative autocontrol may be due to non-specific reactions or may indicate the presence of an alloantibody directed against a high frequency antigen.
 - A positive reaction with all "ID-DiaPanel" test cells and a positive autocontrol may be due to non-specific reactions.
-

- A positive reaction with all "ID-DiaPanel" test cells and the autocontrol but with one or more test cells showing stronger reactions than the autocontrol, may indicate an underlying allo-antibody and further investigation should be undertaken.

IV. Compatibility test

- A negative reaction indicates compatibility of the donor blood with the recipient.
- A positive reaction indicates incompatibility of the donor blood with the recipient, due to presence of antibodies directed against antigens on the donor red cells. Further investigation to identify the antibody specificity should be performed.

REMARKS

As for all procedures covered by GLP rules, the sensitivity of the above procedures should be validated using antibodies of known potency. The Bio-Rad Anti-D reference reagent provides the means to regularly perform controls for all antibody detection procedures (see related package insert).

LIMITATIONS

- a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
- b) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
- c) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
- d) Cells that have become polyagglutinable, due to cryptantigen exposure e.g. T antigen, either *in vivo* or *in vitro* may react with all human sera. Further investigation of such reactions is required.
- e) Bacterial or other contamination of materials used can cause false positive or false negative results.

- f) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
- g) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- h) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
- i) The use of other test cell reagents than "ID-DiaCell" or "ID-DiaPanel" may modify the reaction patterns.
- j) Too heavy or too weak red cell suspensions can cause aberrant results.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performance evaluation studies have been realized with ID-Cards "LISS/Coombs". Results obtained showed that performance characteristics of each application of the card are in conformity with the application field and comparable to other approved products.

Reproducibility

Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards "LISS/Coombs" have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

BIBLIOGRAPHY

1. Coombs, R.R.A., Mourant, A.E. and Race, R.R.: Lancet 1945; 2: 15.
 2. Coombs, R.R.A., Mourant, A.E. and Race, R.R.: Brit. J. Exp. Path 1945; 26: 255.
 3. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
 4. Technical Manual; Roback, J.D. (ed.); 17th ed. 2012; American Association of Blood Banks.
-

PRODUCTS

ID-Card "LISS/Coombs"	4 x 12.	REF 004014
	24 x 12.	REF 004017
	60 x 12.	REF 004016
	112 x 12.	REF 004015

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

DiaMed GmbH
Pra Rond 23
1785 Cressier FR
Schweiz

0123

Changes to the version 05.10 are
shaded gray.

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