**EN INSTRUCTIONS FOR USE****Pot / kit suitable for transport / storage of human tissue****Storage before pack is opened:**

- Store in clean, dry conditions away from direct sunlight.
- Temperature: 10° C to 35° C. RH : 20% to 80%.

**Operational Information:**

- Do not use if package damaged.
- Pots and caps can be stored as low as -80° C for up to five years.
- Frozen plastic is very brittle. Handle with care and avoid shocks or pressure when frozen.
- If to be excessively handled allow to return to room temperature first.
- Do not re-use due to the risk of cross-infection.
- Do not re-sterilise.

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The management system of

# Medfor Products Ltd

Unit 2, Gresham Industrial Estate, Eastern Road,  
Aldershot, Hampshire, GU12 4YD, UK

has been assessed and certified as meeting the requirements of

## Directive 93/42/EEC

on medical devices, Annex V

For the following products

**Sterile and non sterile containers for the transit and storage of human organs and tissues prior to transplantation.**


For placing on the market of Class IIb covered by this certificate, an EC Type Examination Certificate according to Annex III is required.

This certificate is valid from 31 October 2019 until 19 December 2023  
and remains valid subject to satisfactory surveillance audits.

Issue 1. Certified since 17 September 2001.  
and first certified by SGS Belgium on 31 October 2019.

Certification is based on reports numbered GB/PC 200659

Authorised by

Pieter Weterings  
Certification Manager

**SGS Belgium NV, Notified Body 1639**

SGS House Noorderlaan 87 2030 Antwerp Belgium  
t +32 (0)3 545-48-48 f +32 (0)3 545-48-49 [www.sgs.com](http://www.sgs.com)

LPMD5008 - Certificate CE1639 AnnexV\_EN rev. 01

Page 1 of 1

This document is issued by the Company subject to its General Conditions of Certification Services, unless otherwise agreed, accessible at [www.sgs.com/terms\\_and\\_conditions.htm](http://www.sgs.com/terms_and_conditions.htm). Attention is drawn to the limitations of liability, indemnification and jurisdictional issues established therein. The authenticity of this document may be verified at <https://www.sgs.com/en/certified-clients-and-products/certified-client-directory>. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

# Freezable Tissue Storage Kits

Medfor's range of *Freezable Tissue Storage Kits* are intended for the storage of tissue grafts, such as amniotic membrane, bone, cardiovascular, sclera, skin and tendon.

- All caps are tamper-evident
- Each kit consists of a transparent pot inside a transparent pot plus caps
- Each pack consists of a kit double-bagged in peel-back pouches
- Kits can be frozen to -80 ° C for up to five years
- Sterilised by ethylene oxide (EO)
- CE-marked as medical devices for in vivo use
- Medfor is certified to ISO 13485:2016

| Reference  | Volume (ml) |           | Internal Diameter (mm) |           | Height incl. cap (mm) |           | Packaging                           |
|------------|-------------|-----------|------------------------|-----------|-----------------------|-----------|-------------------------------------|
|            | Inner Pot   | Outer Pot | Inner Pot              | Outer Pot | Inner Pot             | Outer Pot |                                     |
| SKA0050Z-K | 50          | 100       | 35                     | 49        | 57                    | 63        | 1 kit per pack<br>10 packs per case |
| SKA0100Z-K | 100         | 300       | 49                     | 72        | 63                    | 83        |                                     |
| SKA0300Z-K | 300         | 500       | 72                     | 87        | 83                    | 110       |                                     |
| SKA0500Z-K | 500         | 1000      | 87                     | 109       | 110                   | 127       | 1 kit per pack<br>8 packs per case  |

Contact us for free evaluation samples and prices

Medfor Products Ltd  
Unit 2, Gresham Industrial Estate  
Eastern Road, Aldershot  
GU12 4YD, UK

[www.medfor.co.uk](http://www.medfor.co.uk) [sales@medfor.co.uk](mailto:sales@medfor.co.uk) T: +44 (0)1252 371181

*This document is for exclusive use of health professionals.*

## DECLARATION OF CONFORMITY

|                             |                                                                                                                                    |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Manufacturer                | Medfor Products Limited<br>Unit 2, Gresham Industrial Estate<br>Eastern Road<br>Aldershot, Hampshire<br>GU12 4YD<br>United Kingdom |
| Representative              | Mr Rob Stephens, Commercial Director                                                                                               |
| Devices                     | Storage Containers and lids - as per attached Device Schedule                                                                      |
| Notified Body               | SGS Belgium NV<br>SGS House Noorderlaan 87<br>2030 Antwerp<br>Belgium                                                              |
| Notified Body Number        | 1639                                                                                                                               |
| EC Certificate number       | GB19/964231                                                                                                                        |
| Classification              | Class IIa                                                                                                                          |
| Classification Rule         | MDD 93/42/EEC Annex IX, Section III, Rule 2                                                                                        |
| Conformity Assessment Route | MDD 93/42/EEC Annex V                                                                                                              |

Medfor Products Ltd. declares, under our sole responsibility, that the devices listed on the attached Device Schedule conform to the relevant provisions of the EC Council Directive 93/42/EEC, dated 14th June 1993 and the relevant amendments of 2007/47/EEC, Annex V.

Medfor Products Ltd agrees to develop, implement and maintain a formally recognised Quality Management System to ensure continued adequacy and efficacy.

Medfor Products Ltd manufactures according to the following harmonised standards:

| Standard                         | Description                                                                                                                                                              |
|----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ISO 13485:2016/EN ISO 13485:2016 | Medical Devices - Quality Management Systems - Requirements for Regulatory Purposes                                                                                      |
| EN ISO 14971:2012                | Medical Devices – Application of Risk Management to Medical Devices                                                                                                      |
| ISO 15223-1:2016                 | Medical Devices - Symbols to be used with Medical Device Labels, Labelling and Information to be supplied – Part 1: General Requirements                                 |
| EN 1041:2008+A1:2013             | Information supplied by the manufacturer of medical devices                                                                                                              |
| ISO 11135:2014                   | Sterilisation of health care products – Ethylene Oxide. Requirements for the development, validation and routine control of a sterilisation process for medical devices. |

|                             |                                                                                                                                                                     |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ISO 11137-1:2015            | Sterilisation of health care products-Radiation-Part 1: Requirements for development, validation and routine control of a sterilisation process for medical devices |
| ISO 11137-2:2015            | Sterilisation of health care products-Radiation-Part 1: Establishing the sterilisation dose                                                                         |
| ISO 11138:2006              | Sterilisation of health care products-Biological Indicators-Part1: General Requirements                                                                             |
| ISO 14644-1:2015            | Cleanrooms and associated controlled environments. Classification of air cleanliness                                                                                |
| ISO 10993-1:2009            | Biological Evaluation of medical devices – Part 1: Evaluation and testing within a Risk Management process                                                          |
| EN ISO 11607-1:2009+A1:2014 | Packaging for terminally sterilised medical devices. Requirements for materials, sterile barrier systems and packaging systems                                      |
| EN ISO 11607-2:2006+A1:2014 | Packaging for terminally sterilised medical devices. Validation requirements for forming, sealing and assembly processes                                            |

Signed .....  
Rob Stephens, Commercial Director

Date: 1 November 2019

## DEVICE SCHEDULE

|              |                                                                             |             |                                                                             |
|--------------|-----------------------------------------------------------------------------|-------------|-----------------------------------------------------------------------------|
| OWE0125D-K   | Pot suitable for freezing to -80° C for transport / storage of human tissue | S248-J      | Pot for transport / storage of human tissue                                 |
| OWE0250D-K   | Pot suitable for freezing to -80° C for transport / storage of human tissue | S253-J      | Pot for transport / storage of human tissue                                 |
| OWE0500D-K   | Pot suitable for freezing to -80° C for transport / storage of human tissue | S296-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| OWE0750D-K   | Pot suitable for freezing to -80° C for transport / storage of human tissue | S348-K      | Pot for transport / storage of human tissue                                 |
| OWE1000D-K   | Pot suitable for freezing to -80° C for transport / storage of human tissue | S400-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| OWE2000D-K   | Pot suitable for freezing to -80° C for transport / storage of human tissue | S401-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| PNA0500E-L   | Pot for transport / storage of human tissue                                 | S402-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| PNC0025D-K   | Pot for transport / storage of human tissue                                 | S403-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| PNC0250D-K   | Pot for transport / storage of human tissue                                 | S406-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| PNC1000D-K   | Pot for transport / storage of human tissue                                 | S411-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| PNF0050D-K   | Pot for transport / storage of human tissue                                 | S412-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| PNF0050D-K2  | Kit for transport / storage of human tissue                                 | S450-K      | Pot suitable for freezing to -80° C for transport / storage of human tissue |
| PNF0050D-K20 | Kit for transport / storage of human tissue                                 | S600-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| PNF0050D-K4  | Kit for transport / storage of human tissue                                 | S610-K      | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| S072         | Kit suitable for freezing to -80° C for transport / storage of human tissue | S631-K      | Accessory for pot for transport / storage of human tissue                   |
| S146-K       | Kit suitable for freezing to -80° C for transport / storage of human tissue | SKA0050Z-K  | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| S147-K       | Kit suitable for freezing to -80° C for transport / storage of human tissue | SKA0100Z-K  | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| S148-K       | Kit suitable for freezing to -80° C for transport / storage of human tissue | SKA0300Z-K  | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| S152-K       | Pot for transport / storage of human tissue                                 | SKA0500Z-K  | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| S212-J       | Kit for transport / storage of human tissue                                 | SKB0050Z-K  | Kit suitable for freezing to -80° C for transport / storage of human tissue |
| S215-K       | Kit for transport / storage of human tissue                                 | SPA0050Z-K5 | Pot suitable for freezing to -80° C for transport / storage of human tissue |
| S216-K       | Kit for transport / storage of human tissue                                 | SPA0100Z-K5 | Pot suitable for freezing to -80° C for transport / storage of human tissue |
| S224-J       | Pot for transport / storage of human tissue                                 |             |                                                                             |
| S227-K       | Kit suitable for freezing to -80° C for transport / storage of human tissue |             |                                                                             |
| S234-J       | Kit for transport / storage of human tissue                                 |             |                                                                             |
| S238-K       | Kit for transport / storage of human tissue                                 |             |                                                                             |
| S241-K       | Kit for transport / storage of human tissue                                 |             |                                                                             |
| S242-K       | Kit for transport / storage of human tissue                                 |             |                                                                             |
| S245-J       | Pot suitable for freezing to -80° C for transport / storage of human tissue |             |                                                                             |

## Freezable Tissue Storage Kits

### Pack Description

Each pack consists of a kit double wrapped in peel-back pouches. Each kit consists of a transparent pot inside a transparent pot plus white caps.

### Application

The transport and storage of tissue grafts, such as amniotic membrane, bone, cardiovascular, sclera, skin and tendon.

### Features

- All caps are tamper-evident
- Kits can be frozen to -80° C for up to five years
- Sterilised by ethylene oxide (EO)
- CE-marked as medical devices for in vivo use

### Orders

This range of products is not available via our website. For pricing, evaluation samples, more information or to place an order, please call us or email [sales@medfor.co.uk](mailto:sales@medfor.co.uk)

## Product Information Leaflet

To download the PDF product information leaflet for this range, go to the [Downloads section](https://www.medfor.co.uk/article/20/brochures-and-other-downloadable-documents/52ab1235014c8c6727e392c3ce8d6e12) (<https://www.medfor.co.uk/article/20/brochures-and-other-downloadable-documents/52ab1235014c8c6727e392c3ce8d6e12>), of our website.

## Regulatory / Quality Assurance

- CE-marked and conforms to Medical Device Directive 93/42/EEC including 2007/47/EC.
- Medfor is approved to ISO 13485:2016.

## Sizes

Available in: 50ml , 100ml , 300ml and 500ml

## Sterility

Available in: Sterile EO

## Labels

Available with: Tissue Bank

## Volume Discounts

This product is eligible for [volume discounts](https://www.medfor.co.uk/article/15/discounts/52ab1235014c8c6727e392c3ce8d6e12) (<https://www.medfor.co.uk/article/15/discounts/52ab1235014c8c6727e392c3ce8d6e12>).

[Login](#)


([https://www.medfor.co.uk/customer?](https://www.medfor.co.uk/customer?xCmd=acshow&jssCart=52ab1235014c8c6727e392c3ce8d6e12)

[xCmd=acshow&jssCart=52ab1235014c8c6727e392c3ce8d6e12](https://www.medfor.co.uk/customer?xCmd=acshow&jssCart=52ab1235014c8c6727e392c3ce8d6e12)).


|                                                                                   |                         |                      |                                                                                                                                                                                                         |                           |                            |                         |                                |                       |  |
|-----------------------------------------------------------------------------------|-------------------------|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------|----------------------------|-------------------------|--------------------------------|-----------------------|--|
|  |                         | TECHNICAL DATA SHEET |                                                                                                                                                                                                         |                           |                            | TDS No 012              |                                |                       |  |
| Issue 06                                                                          |                         |                      |                                                                                                                                                                                                         |                           |                            |                         |                                |                       |  |
| 1                                                                                 | Description             |                      | Freezable Tissue Storage Kit, double bagged, Sterile EO                                                                                                                                                 |                           |                            |                         |                                |                       |  |
| 2                                                                                 | Ref(s)                  |                      | Pot sizes                                                                                                                                                                                               |                           | Blank label on outer pot   |                         | Pre-printed label on both pots |                       |  |
|                                                                                   |                         |                      | 50 ml / 100 ml                                                                                                                                                                                          |                           |                            |                         | SKA0050Z-K                     |                       |  |
|                                                                                   |                         |                      | 100 ml / 300 ml                                                                                                                                                                                         |                           |                            |                         | SKA0100Z-K                     |                       |  |
|                                                                                   |                         |                      | 250 ml/500 ml                                                                                                                                                                                           |                           | S600-K                     |                         |                                |                       |  |
|                                                                                   |                         |                      | 300 ml / 500 ml                                                                                                                                                                                         |                           |                            |                         | SKA0300Z-K                     |                       |  |
|                                                                                   |                         |                      | 500 ml/1000 ml                                                                                                                                                                                          |                           | S610-K                     |                         | SKA0500Z-K                     |                       |  |
| 3                                                                                 | Regulatory              |                      | Conform to the relevant provisions of the EC Council Directive 93/42/EEC, dated 14th June 1993 and all the relevant amendments. Annex V sterility aspects only. Medfor is registered to ISO 13485:2016. |                           |                            |                         |                                |                       |  |
| 4                                                                                 | CE marking              |                      | CE 1639 (Class IIa medical device suitable for storage of human tissue for in vivo use)                                                                                                                 |                           |                            |                         |                                |                       |  |
| 5                                                                                 | Manufacture             |                      | Product assembled in class 7 cleanroom in accordance with BS/EN/ISO 14644-1:2015.                                                                                                                       |                           |                            |                         |                                |                       |  |
| 6                                                                                 | Dimensions and Volume   |                      |                                                                                                                                                                                                         | Height including cap (mm) | Diameter including cap(mm) | Internal neck bore (mm) | Nominal capacity (ml)          | Brimful capacity (ml) |  |
|                                                                                   |                         | SKA0050Z-K           | Inner pot                                                                                                                                                                                               | 57                        | 41.4                       | 34.8                    | 50                             | 50                    |  |
|                                                                                   |                         |                      | Outer pot                                                                                                                                                                                               | 63                        | 55.4                       | 48.6                    | 100                            | 105                   |  |
|                                                                                   |                         | SKA0100Z-K           | Inner pot                                                                                                                                                                                               | 63                        | 55.4                       | 48.6                    | 100                            | 105                   |  |
|                                                                                   |                         |                      | Outer pot                                                                                                                                                                                               | 83                        | 80.2                       | 72.1                    | 300                            | 312                   |  |
|                                                                                   |                         | S600-K               | Inner pot                                                                                                                                                                                               | 87                        | 72.2                       | 64.0                    | 250                            | 257                   |  |
|                                                                                   |                         |                      | Outer pot                                                                                                                                                                                               | 110                       | 95.2                       | 86.5                    | 500                            | 587                   |  |
|                                                                                   |                         | SKA0300Z-K           | Inner pot                                                                                                                                                                                               | 83                        | 80.2                       | 72.1                    | 300                            | 312                   |  |
|                                                                                   |                         |                      | Outer Pot                                                                                                                                                                                               | 110                       | 95.2                       | 86.5                    | 500                            | 587                   |  |
|                                                                                   |                         | S610-K<br>SKA0500Z-K | Inner pot                                                                                                                                                                                               | 110                       | 95.2                       | 86.5                    | 500                            | 587                   |  |
| Outer pot                                                                         | 127                     |                      | 117.8                                                                                                                                                                                                   | 108.7                     | 1000                       | 1096                    |                                |                       |  |
| 7                                                                                 | Material composition    |                      | Jar                                                                                                                                                                                                     |                           | Natural polypropylene      |                         |                                |                       |  |
|                                                                                   |                         |                      | Cap                                                                                                                                                                                                     |                           | White HDPE                 |                         |                                |                       |  |
| 8                                                                                 | Biocompatibility        |                      | Literary review of comparative products and historic usage of this type of product sold by Medfor has established that the kits do not affect the quality or integrity of tissue stored within them.    |                           |                            |                         |                                |                       |  |
| 9                                                                                 | Type of closure         |                      | Tamper-evident screw cap                                                                                                                                                                                |                           |                            |                         |                                |                       |  |
| 10                                                                                | Sterilisation           |                      | Sterilisation has been carried out using ethylene oxide in accordance with ISO 11135-1:2014.                                                                                                            |                           |                            |                         |                                |                       |  |
| 11                                                                                | Labelling               |                      | Case, pack and product labelling in accordance with ISO 15223-1:2016                                                                                                                                    |                           |                            |                         |                                |                       |  |
| 12                                                                                | Temperature range       |                      | -80° C to 75° C. Can be stored to -80° C for up to five years.                                                                                                                                          |                           |                            |                         |                                |                       |  |
| 13                                                                                | Autoclavability         |                      | Not suitable.                                                                                                                                                                                           |                           |                            |                         |                                |                       |  |
| 14                                                                                | Physical resistance     |                      | Flexible to rigid, good mechanical resistance, high chemical resistance                                                                                                                                 |                           |                            |                         |                                |                       |  |
| 15                                                                                | Shelf life (before use) |                      | Packs contents remain sterile for up to 3 years from date of manufacture or until pack is opened or breached, whichever occurs first.                                                                   |                           |                            |                         |                                |                       |  |

|    |                                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |                   |                                                                                                                                                                           |
|----|----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 16 | Pack quantities                        | SKA0050Z-K                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 10 packs per case | Pack description:<br>50 ml pot inside 100 ml pot with both caps separate within inner pack so ethylene oxide unimpeded during sterilisation process. Both pots labelled.  |
|    |                                        | SKA0100Z-K                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 10 packs per case | Pack description:<br>100 ml pot inside 300 ml pot with both caps separate within inner pack so ethylene oxide unimpeded during sterilisation process. Both pots labelled. |
|    |                                        | S600-K                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 10 packs per case | Pack description:<br>250ml pot inside 500ml pot with both caps separate within inner pack so ethylene oxide unimpeded during sterilisation process. Outer pot labelled.   |
|    |                                        | SKA0300Z-K                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 10 packs per case | Pack description:<br>300 ml pot inside 500 ml pot with both caps separate within inner pack so ethylene oxide unimpeded during sterilisation process. Both pots labelled. |
|    |                                        | S610-K                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    | 8 packs per case  | Pack description:<br>500ml pot inside 1000ml pot with both caps separate within inner pack so ethylene oxide unimpeded during sterilisation process. Outer pot labelled.  |
|    |                                        | SKA0500Z-K                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 8 packs per case  | Pack description:<br>500ml pot inside 1000ml pot with both caps separate within inner pack so ethylene oxide unimpeded during sterilisation process. Both pots labelled.  |
| 17 | Storage conditions before packs opened | Temperature: 10°C to 35°C<br>Relative Humidity: 20% to 80%<br>Store in clean, dry conditions away from direct sunlight.                                                                                                                                                                                                                                                                                                                                                                                   |                   |                                                                                                                                                                           |
| 18 | Other information                      | <ul style="list-style-type: none"> <li>Do not use if package damaged.</li> <li>Single Use only. Do not reuse due to risk of cross-infection.</li> <li>Do not irradiate</li> <li>The range of Tissue transport freezing storage pots and Tissue transport storage pots are both DEHP and Latex Free</li> <li>The range does not contain tissues of animal origin nor human blood derivative products.</li> <li>The devices do not require to be connected to other devices in order to be used.</li> </ul> |                   |                                                                                                                                                                           |
| 19 | Product labels                         | SKA0050Z-K – image of label applied to both pots (actual size 60 mm wide x 40 mm high):                                                                                                                                                                                                                                                                                                                                                                                                                   |                   |                                                                                                                                                                           |

SKA0100Z-K – image of label applied to both pots (actual size 60 mm wide x 40 mm high):

S600-K – label size on outer pot is 185 mm wide x 55 mm high (no image)

SKA0300-K – image of label applied to both pots (actual size 90 mm wide x 35 mm high):

S610-K – label size on outer pot is 185 mm wide x 55 mm high (no image).

SKA0500Z-K – image of label applied to both pots (actual size 90 mm wide x 35 mm high):