

hawa

hd 650 DE-V

EcoPak 03/03S

OPERATING INSTRUCTIONS

Valid for device versions of 09/2023

1	INTRODUCTION.....	3
1.1	PREFACE	3
1.2	LEGEND	3
1.3	IMPORTANT NOTICE	4
1.4	CLEANING	4
1.5	SAFETY INSTRUCTIONS	5
2	BEFORE STARTING	7
2.1	INTENDED USE	7
2.2	DESIGN AND FUNCTIONS	8
2.3	INSTALLATION.....	9
3	BASIC FUNCTIONS	10
3.1	SWITCHING ON THE DEVICE	10
3.2	PRINTER SETTINGS	10
3.3	OPERATION	11
3.4	PROCESS VARIABLES	11
3.5	SETTING THE TEMPERATURE CONTROLLER	12
3.6	CONTROL LIGHT FUNCTIONS.....	13
3.7	FUNCTIONS OF THE BUTTONS	14
4	DEVICE SETTINGS.....	15
4.1	SEALING TEMPERATURE ENTRY	15
4.2	SWITCH PRINTER OFF AND ON.....	15
4.3	ENTRY OF A PERSONNEL CODE	16
4.4	DATA ENTRY.....	17
4.5	SELECT PRINT DATA	24
4.6	NUMBER OF PRINTED LABELS	25
4.7	SELECT LABEL SIZES	25
4.8	CONTENTS LABEL INFORMATION FIELD	27
4.9	LABEL INFORMATION – SEALING AND PRINTING	27
5	INTELLIGENTSCAN, CONNECTION OF A BARCODE SCANNER	28
5.1	SEALING SEAM TEST – “SEAL CHECK”	29
6	TROUBLESHOOTING AND MAINTENANCE.....	30
6.1	TROUBLESHOOTING CHECKLIST.....	30
6.2	CUSTOMER SERVICE.....	31
7	ALARM FUNCTIONS AND ERROR DISPLAYS.....	32
7.1	ALARM FUNCTIONS.....	32
7.2	ERROR DISPLAYS	32
7.3	MAINTENANCE/CALIBRATION.....	34
7.4	SPARE PARTS SERVICE	34
7.5	REPLACEMENT PARTS ORDERING – ALLOCATION OF ARTICLE NUMBERS.....	37
7.6	REPLACEMENT PARTS ORDERING - COMPREHENSIVE OVERVIEW.....	38
7.7	REPLACING OF WEARING PARTS	39
7.8	ADJUSTING THE PROCESS PARAMETERS.....	41
7.9	CONTACT PRESSURE	43
7.10	CONFIGURATION OF MOTOR MONITORING.....	45
7.11	ACTIVATE DEFAULT SETTINGS	46
7.12	SETTING THE SERIAL INTERFACE TRANSMISSION RATE (BAUD RATE).....	47
8	TECHNICAL DATA.....	48
8.1	CIRCUIT DIAGRAM AND WIRING DIAGRAM 230V~	48
8.2	CIRCUIT DIAGRAM AND WIRING DIAGRAM 115V~	49
8.3	SPECIFICATIONS	50
9	DECLARATION OF CONFORMITY	51

1 Introduction

1.1 Preface

First of all we would like to thank you for purchasing this sealing device. In these instructions you will find information about using the device, servicing and care.

The sealing device is a microprocessor controlled rotary sealer with a printer for packaging sealable transparent pouches and reels (med. packaging).



Please read these operating instructions carefully before commissioning so that you are familiar with the capabilities of the device and you can make optimum use of its functions.

1.2 Legend

	The exclamation mark in the triangle draws your attention to important notes in the operating instructions, which must absolutely be observed.
	This warning sign refers to measures that could result in danger to human health if they are not observed. It is compulsory to observe it.
	Tips with a hand symbol next to them relate to daily practice.

1.3 Important notice



In accordance with the intended use, the CE marking is displayed based on the following EU directive.

The Medical Device Directive 93/42/EEC and Medical Device Regulation 2017/745 (MDR) is not applicable to sealing devices.

The limit values of IEC 60601-1 must not be applied in repeated electrical inspections.

The manufacturer accepts no liability whatsoever for damage caused by tests in accordance with standards not listed in the declaration of conformity.

In the event of conversion work or interventions to the device undertaken without the express written permission of the manufacturer, the warranty shall be deemed void and any liability for physical or material damage shall be transferred to the operator.

Note

We are constantly improving our products, therefore we reserve the right to modify these operating instructions and the functions described in them.

These operating instructions apply to products from the EcoPak 03/03S series.

1.4 Cleaning

Before cleaning, disconnect the mains plug from the socket and disconnect the device from the power supply with the plug.

Clean the device only with a dry or damp soft cloth and a mild cleaning agent. (E.g.: isopropanol, spirit, etc.) Do not allow any water to find its way into the device.

Caution! Never wet clean the device!

1.5 Safety instructions



1. Our products are in a flawless condition in terms of safety technology when they leave the plant.
2. To maintain this condition, the contents of these safety instructions as well as type plates, labelling and safety instructions attached to the device must be observed while handling the device (transport, storage, installation, commissioning, operation and maintenance).
3. This device is suitable for processing laminated films in the heat-sealing process. See also chapter 2.1 "Intended use".
4. Please check the packaging, and lodge a complaint for any damage with the carrier or parcel service immediately, before installing the device.
5. Before commissioning, ensure that the device does not show any evidence of damage. In case of doubt, contact the manufacturer or a service partner authorised by the manufacturer.
6. Do not operate the device if the power cable or the power plug is damaged. Do not use the device if it does not operate correctly or it is damaged in any way. If the mains cable or the device have been damaged, the device must be repaired by the manufacturer or by one of the manufacturer's authorised service partners.
7. The device must be connected using the mains cable included in the scope of delivery to a protective contact socket with a stable voltage. Operation on IT networks is not permitted.
8. Place the device on a stable base.
9. The device must not be installed or operated in potentially explosive areas.
10. If the sealing device is brought directly from a cold environment into a warm environment, condensation may form. Wait until temperature equalisation has taken place.
Starting up the device when it contains condensation causes danger to life!
11. Repairs and the replacement of wear parts / spare parts must be performed only by the manufacturer or by one of the manufacturer's authorised service partners.
12. Switch off the device when it is not in use, or remove the power plug from the socket.
13. Before cleaning: Disconnect from the mains! Clean the device only with a dry or damp soft cloth and a mild cleaning agent. Do not allow any water to find its way into the device. Caution! Never wet clean the device!
14. Do not insert pointed or flat items into the import slot of the device. This can result in damage to the device and instruments.
15. Do not insert items into the louvres of the device. You may receive an electric shock or the device could be damaged.
16. Do not use the device if you have any doubts about its safety.
17. The device must not be installed or operated by persons under 16 years of age.

18. The device must not be operated unsupervised.
19. It is forbidden to operate the device under the influence of drugs or alcohol.
20. Keep hair, clothing and gloves away from moving parts! Loose clothing, jewellery or long hair can be caught by moving parts.



21. The device contains valuable materials that can be recycled and reused. The device should therefore be disposed of at a public disposal facility near you. The device has been labelled in accordance with Directive 2012/19/EU (WEEE) on waste electrical and electronic equipment. This directive governs the return and recycling of scrap equipment within the EU.

2 Before starting

2.1 Intended use

SEALING MATERIALS

Sealable paper pouches in accordance with EN ISO 11607-1/EN 868-4	x
Sealable pouches and tubes in accordance with EN ISO 11607-1/EN 868-5 made of film and paper as per EN 868-3	x
Sealable pouches and tubes in accordance with ISO EN 11607-1/EN 868-5 made of film and uncoated materials made of polyolefins as per EN 868-9 (e.g. Tyvek ^{®1})	x
Sealable pouches and tubes in accordance with ISO 11606-1/EN 868-5 made of PP fleece or PP non-woven	x ²


NON-SEALABLE MATERIALS

Soft PVC films
Polyamide films
Coated HDPE
Aluminium-laminate film
Polyethylene films
Hard PVC films
Polypropylene films

¹ Tyvek[®] is a registered trademark of E.I. du Pont Nemours.

² Release and / or test required

2.2 Design and functions



Sealing process sequence

- Step 1:** After the med. packaging has been inserted, the feed is automatically switched on.
- Step 2:** The med. packaging is now fed and the sealing seam area is heated up to the set sealing temperature by the heating units located at the top and bottom.
- Step 3:** The sealing seam, which is now heated, is pressed together by the sealing rollers and sealed.
- Step 4:** The finished med. packaging is transported to the extraction side.
- Step 5:** If no item to be sealed is fed in, the feed switches off after approximately 30 seconds.

2.3 Installation



The device must not be installed or operated in potentially explosive areas.


Only use sockets that are equipped with a protective conductor and where the mains voltage is stable.



The device may only be installed in a dry environment. Heavy dust, steam, dripping water or splashing water impair the function of the device.


Please make sure that the operating voltage corresponds to the specifications on the device's type plate.

- 1 Please do **not** transport the device on the peel edge setting.
- 2 The distance from the device to a wall must be at least 200mm!



3 Basic functions

3.1 Switching on the device



Plug the mains cable into the mains connection.

Switch on the device with the mains switch in position "1".

Set the required sealing temperature on the temperature controller as described in Chapter 3.3.


The device is heated up and ready for operation as soon as the set sealing temperature is displayed.

3.2 Printer Settings

All printer functions are only available in conjunction with a corresponding printer
ValiPrint + print server possible.
Please feel free to contact us about this!



3.3 Operation



- Schritt 1: Set the required sealing edge width. After loosening the knurled screw, the lower insertion plate can be continuously adjusted for sealing edge widths of 0 - 30mm.
- Schritt 2: Insert packaging into the device from the left via the insertion plate.
- Schritt 3: Remove the sealed package on the outlet side and allow to cool down briefly.



Check of the sealing seam

If leaks appear, the sealing temperature must be increased. If the film melts, the set temperature is too high.

As per DIN 58953 -7, the suitable sealing temperature is to be determined by test sealing.

3.4 Process variables

Sealing temperature

The temperature is monitored electronically by means of a temperature sensor. If this varies from the set value by 5°C (requirements set out in DIN 58953-7), the drive is locked.






3.6 Control light functions



Function				
Sealing device switched on Heating-up phase	flashes	off	off	off
Sealing temperature = ± 5°	on	off	on	off
Sealing temperature <> ± 5°	off	on	off	on
Standby	off	on	off	on
Heating up after standby or changing the set temperature	flashes	off	off	on








3.7 Functions of the buttons



			
<p>Menu level 1</p> <p>Activation of menu level 2 Press 3 s</p> <p>Activation of menu level 3 Press 7 s</p>	<p>Switch printer off/on</p> 	<p>Enter personal identification</p>	<p>Activation of "seal check"</p>
<p>Menu level 2</p> <p>Entering sealing temperature</p>	<p>Temperature value + 1</p>	<p>Temperature value + -1</p>	<p>Confirm entry</p>
<p>Menu level 3</p> <p>3.1 Sealing parameter view</p> <p>3.2 Print data configuration</p> <p>3.3 Data entry</p>	<p>Changeover 3.1 - 3.2 - 3.3</p> <p>on</p> <p>Entry value +1</p>	<p>Changeover 3.1 - 3.2 - 3.3</p> <p>off</p> <p>Entry value -1</p>	<p>Confirm entry</p>








4 Device settings

4.1 Sealing temperature entry


<p> S</p> <p>Press button 3s</p> <p>Change set temperature</p> <p> </p>	<p>Display for 1s</p>  <p>Display the current temperature setpoint</p>  	<p> Confirm entry</p>
---	---	--

4.2 Switch printer off and on (Just with external printer and printserver)







































<p> ▲</p> <p>Press the button</p> <p>Switch printer off and on</p> <p> </p>	<p>Display for 1s</p>  <p>Display the current setting</p> <p>Printer switched on</p>  <p>Printer switched off</p> 	<p> Confirm entry</p>
--	---	--

4.3 Entry of a personnel code












4.4 Data entry






<p>Activate</p>  <p>Press button 7s</p> <p>Toggle to data entry</p>   <p>Key lock</p> <p>Activated Enter 1 - 9999</p> <p>Deactivated Enter 0</p>  	   	 Confirm selection <p>Select data</p>    Confirm selection  Confirm entry
<p>Date and time</p> <p>Enter the day</p>   <p>Enter month</p>   <p>Enter year</p>   <p>Enter hour</p>   <p>Enter minute</p>  	     	 Confirm selection  Confirm entry  Confirm entry  Confirm entry  Confirm entry  Confirm entry







<p>Expiry date</p> <p>Enter the day</p> <p>Enter month</p> <p>Enter year</p>		Confirm selection Confirm entry Confirm entry Confirm entry
--	--	--








<p>Batch Enter 0000 - 9999</p>		Confirm selection <p> On this display, an alphanumeric code has been entered via the barcode scanner</p> Confirm entry
---	--	---

<p>Package content quantity Enter 0 - 99</p>		Confirm selection Confirm entry
---	--	------------------------------------







<p>Batch counter Enter 0 - 9999</p>  <p>Select counting direction</p> <p>If the counting direction is selected downwards, a signal sounds on reaching the value 0</p> <p>Switch off batch counter</p>	  <p>Counting direction upwards</p>  <p>Counting direction downwards</p>  <p>Switch off batch counter</p> 	 Confirm selection  Confirm entry  Confirm selection
---	---	---

<p>Standby Enter 0 – 60min</p>  <p>If the motor does not start until the set time is reached, the temperature is automatically set to 50°C / 122°F.</p>	  <p>By briefly inserting a package, the temperature returns to it previously set value. Once this is reached, the sealing device is ready for use again</p>	 Confirm selection  Confirm entry
---	--	--


<p>Operating data View operating hours and the absolute batch counter</p> <p>Select view</p> 	 <p>Operating hours</p>  <p>Absolute batch counter</p> 	 Confirm selection  Confirm selection
---	--	--

<p>Measurement units</p> <p>Select measurement units</p> <div style="display: flex; justify-content: space-around;">   </div> <p>Europe Temperature in °C</p> <p>USA Temperature in °F</p>	<div style="text-align: center;">  <p>Measurement units - Europe</p>  <p>Measurement units - USA</p>  </div>	<div style="text-align: center;">  <p>Confirm selection</p>  <p>Confirm selection</p> </div>
---	---	---



<p>Character width</p> <div style="display: flex; justify-content: space-around;">   </div> <p>Enter 0 - 2</p> <p>Depending on the selected setting, the characters are printed in different widths</p> <p>Enter A Automatic adjustment of the character size according to the width of the packaging and the length of the print line.</p> <div style="margin-top: 20px;"> <div style="border: 1px solid black; padding: 2px; width: 100px; text-align: center;">PNR</div> <div style="border: 1px solid black; padding: 2px; width: 100px; text-align: center; margin-top: 5px;">PNR</div> </div>	<div style="text-align: center;">   </div>	<div style="text-align: center;">  <p>Confirm selection</p>  <p>Confirm entry</p> </div>
---	--	---

**Monitoring
Personnel number**
Enter 0 – 60min






When entering a time >0, personal number monitoring is activated and the motor can only ever be started if the personal number value is 1-9999


If the set time is reached, the personal number is automatically set to 0

If the personal number is 0 and a package is inserted, a prompt appears







By entering a personal number 1-9999, the lock on the drive motor is released and the prompt disappears









Confirm selection



Confirm entry

<p>Date format</p> <p>Selection of the date format</p> <p>   </p> <p>1 DD.MM.YYYY</p> <p>2 MM.DD.YYYY</p> <p>3 YYYY.MM.DD</p> <p>4 YYYY.MM</p> <p>5 DD.MM.YY</p> <p>6 MM.DD.YY</p> <p>7 YY.MM.DD</p>	 	<p> Confirm selection</p> <p> Confirm selection</p>
---	--	--

<p>Time format</p> <p>Select time format</p> <p>   </p> <p>24 13:26</p> <p>12 01:26 PM</p>	<p>24-hour time format</p>  <p>12-hour time format</p> 	<p> Confirm selection</p> <p> Confirm selection</p>
---	--	--




4.5 Select print data

Activate




Press button 7s

Switching of print data




The printing sequence is set and cannot be changed. Activated print data are always printed in this order

- Sterilisation date**
- Time**
- Expiry date**
- Batch number**
- Personnel identification**
- Batch counter**
- Sterilisation type**
- Package content quantity**
- Text**
(via barcode list with connected barcode scanner)




Sterilisation date



Time



Expiry date



Batch number




Personnel identification



Batch counter






Sterilisation type









Confirm selection

Select data to be printed









4.6 Number of printed labels




<p>Enter 1-10 Factory setting 1</p> <p>Select quantity</p> <p> </p>	 	<p> Confirm selection</p> <p> Confirm entry</p>
---	--	--

4.7 Select label sizes




<p>Enter 1-5 Factory setting 1</p> <p>Select sizes</p> <p> </p> <p>Size 1 55x31 mm</p> <p>Size 2 60x44 mm</p> <p>Also see section 3.5.3</p>	 	<p> Confirm selection</p> <p> Confirm entry</p>
---	---	---



Example

Printing the sterilisation date and the expiration date



5 IntelligentScan, connection of a barcode scanner

The following entries and functions can be implemented using a hm 980 BR barcode scanner (item number 1.421.018) connected to the "Intelligent Scan" interface (see page 9), and relevant barcode lists:

Entries and functions using the controller or the barcode scanner

Entries

- Sealing temperature entry
- Personnel code entry
- Batch number entry
- Quantity entry of the packaging contents
- Sterilisation type selection
- Default setting of the batch counter
- Character width selection
- Print data selection

Functions

- Switching the printer off or on
- Activating/deactivating the standby function
- Monitoring the personal number
- Switching off or on
- Seal check activation

Entries and functions only using the barcode scanner

Entries

- Entry of a 10-digit alphanumeric personnel code
- Entry of a 10-digit alphanumeric batch designation
- Entry of an alphanumeric text
- Expiry dates in 1,3,6,9,12,24 and 60 months

Functions

- Switching the item counter off or on



The hm 980 BR barcode scanner (item number 1.421.018) is supplied with a CD (item number 1.490.028) enclosed, facilitating the generation and recording of the barcode lists on a PC.



Please only use barcode scanners approved by hawo. hawo accepts no liability for damage caused by connecting and using other barcode scanners.

For further information, please contact your authorised service partner or the hawo Service Hotline: +49 (0) 6261 9770 0

5.1 Sealing seam test – “Seal Check”

Test the critical process parameters (temperature, contact pressure and sealing time) with "SEAL CHECK".





Seal Check seal indicators are not suitable for packaging with side folds.

This test should be performed before and after the daily working process and/or before/after each batch and can be documented by routine filing of the print-out.


Additional use of the SEAL CHECK sealing indicator in combination with the SEAL CHECK function of the sealing device is recommended.

Before the test, the sealing device must be ready for use and the set temperature must have been reached.

<p>Activate seal check</p> <p>Place the seal indicator into the package</p> <p>Carry out seal procedure</p>	<p>Press and hold the key for 2 s until the Seal Check function is activated.</p>  <p>Then release the key</p>	
--	--	--

6 Troubleshooting and maintenance

6.1 Troubleshooting checklist

 The troubleshooting suggestions marked with a * may be carried out only by the manufacturer or an authorised service partner appointed by the manufacturer.

Malfunction	Possible cause	Remedy
Machine fails to switch on No data on the display	Power supply -Power cable not plugged in	Check power supply If necessary, plug into a different power outlet
	-Power cable defective	Replace power cable
	Fuse	Replace the fuse* ! If the fuse blows again, it is imperative to have the machine tested.
	Temperature controller faulty	Replace temperature controller*
Machine fails to heat up	Set temperature is too low	Increase set temperature (see 3.3.1)
	Temperature limitation activated	Switch off the machine and allow it to cool down ! If it still trips it is imperative to have the machine tested
	Temperature sensor	Replace temperature sensor*
	Heating cartridge	Check heating cartridges and replace if necessary*
	Temperature controller faulty	Replace temperature controller*
No transport	Transport belt -damaged -no transport	Replace transport belt Check belt tension
	Motor sensor	Replace light barrier*
	Motor Temperature controller faulty	Replace motor*
	Temperature controller faulty	Replace temperature controller*

Malfunction	Possible cause	Remedy
Uneven material feed or loud running noise	Transport belt guide	Renew PTFE belt on guide rail (see pg. 36)
	Transport belt -damaged -no transport	Replace transport belt Check belt tension
	Motor	Replace motor*
Sealing seam will not hold	Temperature too low	Increase temperature
	Contact pressure too low	Readjust the contact pressure of the sealing roller or replace the sealing roller*
	Sealing die Distance between the sealing dies too great	Set the distance of the sealing dies to 0.5 mm*
Sealed seam distorted	Pressure applied too high	Readjust the contact pressure of the sealing roller or replace the sealing roller*
Paper side of the package is discoloured or side fold shrunk	Temperature too high	Reduce to temperature (see pg. 14 3.4.1)
No printing or printing incomplete	Ink ribbon	Ink ribbon not inserted properly Replace ink ribbon. (see pg. 35)
	Print head	Replace print head*
	Printer control faulty	Replace printer control*
	Switching power supply faulty	Replace switching power supply
Printing too faint	Ink ribbon	Replace ink ribbon.
	Print head	Readjust print head*
	Paper pressure roller	Adjust paper pressure roller*

6.2 Customer service




Your hawa customer service team is available from Mon - Fri 8 a.m. – 4:30 p.m. and on +49 (0)6261-9770-0.


You are also welcome to send questions to the following e-mail address: service@hawa.com


7 Alarm functions and error displays


7.1 Alarm functions


<p>Batch counter, set with counting direction downwards, has reached the value 0</p>	
<p>Cancelling the alarm</p>	
<p>Set batch counter to a value > 0 or set the counting direction of the batch counter to up See page 19</p>	

7.2 Error displays

 The troubleshooting suggestions marked with a * may be carried out only by the manufacturer or an authorised service partner appointed by the manufacturer.

<p>Sealing temperature outside Tolerance</p>	
<p>Possible cause</p>	
<p>Temperature sensor defective</p>	<p>Replace temperature sensor*</p>
<p>Controller defective</p>	<p>Replace controller*</p>
<p>SST module faulty</p>	<p>Replace SST module*</p>

<p>Contact pressure out of tolerance</p>	
<p>Possible cause</p>	
<p>DMS module not calibrated</p>	<p>Recalibrate DMS module*</p>
<p>DMS module faulty</p>	<p>Replace DMS module*</p>

<p>Sealing speed outside tolerance</p>	 <p>The diagram shows a control panel with a red LED display showing 'E In E'. Above the display are three indicator lights labeled 'okay', 'stby', and 'Parameter'. The 'Parameter' light is red and labeled 'lights up'. Below the display are four buttons labeled 'S', '▲', '▼', and '←'. An arrow labeled 'flashes' points to the display, and another arrow labeled 'lights up' points to the 'Parameter' light.</p>	
	<p>Possible cause</p>	<p>Remedy</p>
	<p>Motor defective</p>	<p>Replace motor*</p>

7.3 Maintenance/calibration



Like all technical devices, your device is also subject to technical wear. In order to guarantee continuous operational readiness, your device should be inspected regularly by a competent person and serviced and calibrated at least once per year by the manufacturer or by one of the manufacturer's authorised service partners.

Maintenance cycle	Ink ribbon	PTFE tape Guiding die	Pressure roller	Toothed belt	Distance of Sealing die	Calibration of process variables
At least every 3 months						
Depending on usage, at least once a year						

Legend:



Check



Replace



Adjust



Measure

7.4 Spare parts service



Simply order parts by fax:

- Please copy the following pages according to the part required
Page 31: Maintenance and wear parts
Page 32: Spare parts
- Enter the device number. _____ →
- Enter the device type. _____ →
- Enter address, fax number and order number.
- Mark items required.
- Enter quantity required.
- Sign order.
- Fax order.

S/N: 555555
Type: hd 650 DE-V

Bitte diese Daten bei Ersatzteilbestellung angeben.
Please state this data when ordering spares.
Veuillez indiquer ces données en cas de commandes
des pièces de rechange.

hawa Gerätebau GmbH
D-74847 Obriegheim
Tel. 06261-62016
Fax. 06261-62015
Made in Germany

2 000 004

To:

Sender:

Fax no.

Your order no. _____		Date _____	
Device type _____		Serial number _____	
<input checked="" type="checkbox"/>	Designation	Art. No.:	Qty.
<input type="checkbox"/>	Ink ribbon, black	6.813.104	
<input type="checkbox"/>			
<input type="checkbox"/>	PTFE strip on upper guide rail	6.105.178	
<input type="checkbox"/>	PTFE strip on lower guide rail	6.105.177	
<input type="checkbox"/>	PTFE strip heating die	6.105.125	
<input type="checkbox"/>	Plastic pressure roller	2.230.008	
<input type="checkbox"/>	Toothed belt, drive	6.271.002	
<input type="checkbox"/>	Toothed belt, transport sealing material	6.271.001	
<input type="checkbox"/>	Heating cartridge	6.536.024	
<input type="checkbox"/>	Upper sealing die assembly	1.616.049	
<input type="checkbox"/>	Complete lower sealing die	1.616.050	
<input type="checkbox"/>			
<input type="checkbox"/>			

To:


Sender:

Fax no.


Your order no. _____		Date _____	
Device type _____		Serial number _____	
<input checked="" type="checkbox"/>	Designation	Art. No.:	Qty.
<input type="checkbox"/>	Temperature controller 100 - 240V	6.564.044	
<input type="checkbox"/>	transformer	6.531.015	
<input type="checkbox"/>	SST module	1.461.014	
<input type="checkbox"/>			
<input type="checkbox"/>			
<input type="checkbox"/>	Motor optical sensor	1.561.010	
<input type="checkbox"/>	Gear motor 230V	1.212.023	
<input type="checkbox"/>	Lüfter 230V	6.212.029	
<input type="checkbox"/>	Lüfter 115V	6.212.030	
<input type="checkbox"/>			
<input type="checkbox"/>	Reset the temperature limiter	6.564.018	
<input type="checkbox"/>	Temperature sensor	6.564.023	
<input type="checkbox"/>			

7.5 Replacement parts ordering – allocation of article numbers


Heating die





Upper guide rail



Lower guide rail



7.6 Replacement parts ordering - comprehensive overview



7.7 Replacing of wearing parts


Maintenance information

! Please use **only genuine replacement parts**


Replacement of PTFE tape on guide rail

→ Switch off the sealing device and DISCONNECT THE POWER PLUG!

- Opening the housing
- Remove mounting screws **1** for the upper guide rail and remove the guide rail **or**
- Remove mounting screws **2** for the lower guide rail and remove guide rail
- Remove mounting screws **3** and detach the PTFE strip
- Pull backing foil off of new PTFE strip and glue new PTFE strip on straight and without wrinkles
- Fasten PTFE strip with screw **3**
- Install the guide rails.

 When installing the upper guide rail, push the die down before fastening so that the gap between the screw head and the rail is 1mm on both sides. This ensures the correct contact pressure for the guide rail.

- Close housing




Maintenance information

Replacing the pressure roller

➔ **Switch off the sealing device and DISCONNECT THE POWER PLUG!**


- Opening the housing
- Remove mounting screws **1** for the upper guide rail and remove the guide rail.
- Unscrew pressure adjustment screw **2** approx. 5 mm
- Loosen mounting screw **3** and pull the pressure roller completely out of the holder
- Detach snap ring **4** and remove the pressure roller
- Install the new pressure roller and fasten with snap ring **4**
- Place the pressure roller fully in the holder, aligning it centrally with the lower roller
- Tighten the mounting screw **3**
- Adjust contact pressure by screwing in adjustment screw **2** according to calibration instructions on page 41




-
- Install guide rail



When installing the upper guide rail before fastening, push the rail as far down as possible, before final fixing with the mounting screw, **1** so that the gap between the screw head and the rail is 1mm on both sides. This ensures the correct contact pressure for the guide rail.




Distance between screw head and rail = 1mm



- Close housing

7.8 Adjusting the process parameters

 After adjusting, the sealing device must remain switched on for another 10 seconds!

Temperature control



Adjustment of temperature control must always be performed after replacing a heating cartridge, after replacing the temperature sensor and after replacing the temperature controller.

The temperatures **120°C and 200°C** are measured one after another and the difference between the set value and the actual value is corrected.


After a set temperature has been reached it is stabilised for 120s. After the 120 seconds have elapsed, the measured temperature value is entered using the temperature measuring instrument

Procedure


Insert the temperature measuring instrument's sensor between the sealing dies from the left infeed side



Press both buttons at the same time and switch on the sealing device




Power switch











7.9 Contact pressure

Press both buttons at the same time and switch on the machine

Adjustment of contact pressure must always be performed after replacing the temperature controller, after replacing the DMA module and after replacing the pressure roller.




<p>Select contact pressure setting</p> <div style="display: flex; justify-content: space-around;">   </div>	 	<div style="display: flex; align-items: center; justify-content: center;">  <div style="margin-left: 10px;">Confirm selection</div> </div>
--	---	---












<p>Procedure</p> <p>Turn screw S to the left until the displayed value no longer changes</p> <p>Using the potentiometer P1, set the 0-point of the DMS module to the value 005</p> <p>By turning the screw S to the right, set the contact pressure to 100N</p>	 <p>The diagram shows a cross-section of the device. A potentiometer labeled 'P1' is located on the top surface. A screw labeled 'S' is located on the front panel, pointing to a specific adjustment point.</p>	 <p>The first display shows '0000'. The second display shows '0005'. The third display shows '1000'.</p> <p> Accept setting</p>
---	---	---

7.10 Configuration of motor monitoring

Press both buttons at the same time and switch on the machine


Configuration of the motor monitoring must always be performed after replacing the temperature controller, after replacing the motor and after replacing the tacho generator



<p>Select configuration of motor monitoring</p> <p> </p> <p>Select frequency of motor voltage</p> <p> </p> <p>The motor starts up and the measured sealing speed is displayed for 5 seconds</p>	<p></p> <p> Power switch</p> <p>50Hz </p> <p>60Hz </p> <p></p>	<p> Confirm selection</p> <p> Confirm selection</p>
---	--	--


7.11 Activate default settings

Press both buttons at the same time and switch on the machine




	Expiry date	6 months	
	Personal identification	Empty	is not printed
STERILE	Sterilisation type	STEAM	
Σ	Package content quantity	1	
LOT	Batch identification	Year and batch counter	
	Text	is not printed	


Select default setting activation



Display flashes



Confirm selection



Confirm activation

7.12 Setting the serial interface transmission rate (baud rate)


Press both buttons at the same time and switch on the sealing device

Power switch

<p>Select baud rate setting</p> <div style="display: flex; justify-content: space-around;"> <div style="border: 1px solid black; border-radius: 50%; width: 40px; height: 40px; display: flex; align-items: center; justify-content: center;">▲</div> <div style="border: 1px solid black; border-radius: 50%; width: 40px; height: 40px; display: flex; align-items: center; justify-content: center;">▼</div> </div> <p>Select baud rate</p> <p>1200Bd – 57600Bd</p> <div style="display: flex; justify-content: space-around;"> <div style="border: 1px solid black; border-radius: 50%; width: 40px; height: 40px; display: flex; align-items: center; justify-content: center;">▲</div> <div style="border: 1px solid black; border-radius: 50%; width: 40px; height: 40px; display: flex; align-items: center; justify-content: center;">▼</div> </div>	<div style="border: 1px solid black; width: 100%; height: 30px; display: flex; align-items: center; justify-content: center; margin-bottom: 10px;"> 8.5 </div> <div style="border: 1px solid black; width: 100%; height: 30px; display: flex; align-items: center; justify-content: center; margin-bottom: 10px;"> BAUD </div> <p style="text-align: center;">Example: Baud rate 57600Bd</p> <div style="border: 1px solid black; width: 100%; height: 30px; display: flex; align-items: center; justify-content: center;"> 576 </div>	<div style="border: 1px solid black; border-radius: 50%; width: 60px; height: 60px; display: flex; align-items: center; justify-content: center; margin-bottom: 20px;"> ← </div> <p style="text-align: center;">Confirm selection</p> <div style="border: 1px solid black; border-radius: 50%; width: 60px; height: 60px; display: flex; align-items: center; justify-content: center;"> ← </div> <p style="text-align: center;">Confirm selection</p>
---	---	--


8 Technical data

8.1 Circuit diagram and wiring diagram 230V~



0	DMS – module	1.410.018
1	opto-sensor	1.561.010
2	gear motor 230V	1.212.023
3		
4	heating cartridge	6.564.024
5	over temperature limit switch	6.564.018
6	SST module	1.461.014
7	transformer 115V/230	6.531.015
8		
9		
10		
1 1		
1 2	temperature controller	6.564.044
1 3	temperature sensor	6.564.023
1 4		

8.2 Circuit diagram and wiring diagram 115V~



0	DMS – module	1.410.018
1	opto-sensor	1.561.010
2	gear motor 230V	1.212.023
3		
4	heating cartridge	6.564.024
5	over temperature limit switch	6.564.018
6	SST module	1.461.014
7	transformer 115V/230	6.531.015
8		
9		
10		
11		
12	temperature controller	6.564.044
13	temperature sensor	6.564.023
14		

8.3 Specifications

Connection data

Mains connection	[V]	115 / 230
Mains frequency	[Hz]	50 / 60
power consumption	max. heating up average [W]	500 390 ~ 160
Mains fuse 115V / 230V	[A]	5A / 2A

Mechanical system

Dimensions	Length	[mm]	465
Including Infeed section	Width		255
	Height		155
Housing			Metal, powder-coated
Weight	230V~ 115V~	[kg]	14,7 15,7
Seal distance from edge		[mm]	0 – 35
Sealing seam width		[mm]	12 +/- 1
Sealing system			Seal Peak
Sealing seam length		[mm]	Unlimited
Distance from medical product		[mm]	>30 (as per DIN 58953-7)

Process variables/sealing parameters

Sealing temperature	max.	[°C]	220
Tolerance for sealing temperature		[°C]	±5
Throughput speed [fixed]		[m / min]	10
Temperature standard tolerance		[%]	±2
contact pressure		[N]	100 (factory setting)
Contact pressure switch-off tolerance		[%]	±20



Electronics and communication systems

System		electronic
interface	RS232	yes
Electrical protection class		1

Environmental parameters

Ambient temperature		[°C]	5-25
Heat output		[kJ/s]	0.1
Noise intensity acc. to Machinery Directive 2006/42/EC Appendix I 1.7.4.2 u.)		[dB/ A]	<70
Relative humidity		[%]	30-80 non-condensing

9 Declaration of conformity

 74847 Obrigheim / Germany	Konformitätserklärung – Declaration of Conformity Déclaration "CE" de Conformité Declaración de conformidad de la C.E. Dichiarazione di conformità - Declaração de conformidade		9.694.048C																																																																																											
	Gültig ab: 24.07.2023 Valid from: 24.07.2023		Seite 1/1 Version 1.01																																																																																											
<p>Hiermit erklären wir, daß die Folienschweissmaschinen: Herewith we declare that the Foil sealing unit: Par la présente, nous déclarons que la gamme de Soudeuse de films plastique: Por la presente certificamos que las máquinas embolsadoras modelos: Dichiariamo con la presente che le macchine per saldatura di fogli: Por este meio se declara que as máquinas de selagem de folhas de plástico:</p> <p style="text-align: center;">hd 650 DE-V ECOPAK 03/03S</p> <p>folgenden einschlägigen Bestimmungen und harmonisierten Normen entsprechen: complies with the requirements of the following regulations and harmonised standards: corresponde aux dispositions suivantes et standards harmonisés: objeto de esta Declaración cumple con las siguientes disposiciones: Sono conformi alle seguenti disposizioni in materia nonché alle seguenti norme armonizzate: correspondem às seguintes determinações e normas harmonizadas:</p> <table border="0"> <tr> <td>EG - Maschinenrichtlinie</td> <td></td> <td>2006/42/EG</td> <td></td> <td></td> </tr> <tr> <td>Machinery directive</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Directive "CE" rel. aux machines</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Directiva de Maquinaria de la CE</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Direttiva CE sulle macchine nella versione</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Directiva da UE relativa a maquinaria</td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>EMV-Richtlinie</td> <td>EMC-directive</td> <td>2014/30/EU</td> <td></td> <td></td> </tr> <tr> <td>Directive CEM</td> <td>Directiva de CEM</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Diretiva CEM</td> <td>Directiva CEM</td> <td></td> <td></td> <td></td> </tr> <tr> <td>WEEE-Richtlinie</td> <td>WEEE-directive</td> <td>2012/19/EU</td> <td></td> <td></td> </tr> <tr> <td>Directive WEEE</td> <td>Directiva de WEEE</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Diretiva WEEE</td> <td>Directiva WEEE</td> <td></td> <td></td> <td></td> </tr> <tr> <td>RoHS-Richtlinie</td> <td>RoHS-directive</td> <td>2015/863/EU</td> <td></td> <td></td> </tr> <tr> <td>Directive RoHS</td> <td>Directiva de RoHS</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Diretiva RoHS</td> <td>Directiva RoHS</td> <td></td> <td></td> <td></td> </tr> <tr> <td>Harmonisierte Normen</td> <td>Harmonized standards</td> <td>EN ISO 12100:2010</td> <td>EN 60204-1:2018</td> <td></td> </tr> <tr> <td>Standard harmonise</td> <td>Las normas armonizadas</td> <td>EN ISO 13857:2019</td> <td>EN IEC 61000-6-1:2019</td> <td></td> </tr> <tr> <td>Norme armonizzate</td> <td>Normas harmonizadas</td> <td>EN IEC 63000:2018</td> <td>EN IEC 61000-6-3:2021</td> <td></td> </tr> </table> <p>Verantwortliche Person für die Technischen Unterlagen siehe unten Responsible person for technical documentation see below La personne responsable pour la documentation technique est mentionnée au-dessous</p> <div style="text-align: center;">  hawo GmbH 74847 Obrigheim, Germany T +49 (0) 6261 9770-0 </div> <p>Torsten Ehrhardt info@hawo.com Prokurist / authorized officer</p> <p>hawo GmbH, Obere Au 2, D-74847 Obrigheim, Germany</p>					EG - Maschinenrichtlinie		2006/42/EG			Machinery directive					Directive "CE" rel. aux machines					Directiva de Maquinaria de la CE					Direttiva CE sulle macchine nella versione					Directiva da UE relativa a maquinaria					EMV-Richtlinie	EMC-directive	2014/30/EU			Directive CEM	Directiva de CEM				Diretiva CEM	Directiva CEM				WEEE-Richtlinie	WEEE-directive	2012/19/EU			Directive WEEE	Directiva de WEEE				Diretiva WEEE	Directiva WEEE				RoHS-Richtlinie	RoHS-directive	2015/863/EU			Directive RoHS	Directiva de RoHS				Diretiva RoHS	Directiva RoHS				Harmonisierte Normen	Harmonized standards	EN ISO 12100:2010	EN 60204-1:2018		Standard harmonise	Las normas armonizadas	EN ISO 13857:2019	EN IEC 61000-6-1:2019		Norme armonizzate	Normas harmonizadas	EN IEC 63000:2018	EN IEC 61000-6-3:2021	
EG - Maschinenrichtlinie		2006/42/EG																																																																																												
Machinery directive																																																																																														
Directive "CE" rel. aux machines																																																																																														
Directiva de Maquinaria de la CE																																																																																														
Direttiva CE sulle macchine nella versione																																																																																														
Directiva da UE relativa a maquinaria																																																																																														
EMV-Richtlinie	EMC-directive	2014/30/EU																																																																																												
Directive CEM	Directiva de CEM																																																																																													
Diretiva CEM	Directiva CEM																																																																																													
WEEE-Richtlinie	WEEE-directive	2012/19/EU																																																																																												
Directive WEEE	Directiva de WEEE																																																																																													
Diretiva WEEE	Directiva WEEE																																																																																													
RoHS-Richtlinie	RoHS-directive	2015/863/EU																																																																																												
Directive RoHS	Directiva de RoHS																																																																																													
Diretiva RoHS	Directiva RoHS																																																																																													
Harmonisierte Normen	Harmonized standards	EN ISO 12100:2010	EN 60204-1:2018																																																																																											
Standard harmonise	Las normas armonizadas	EN ISO 13857:2019	EN IEC 61000-6-1:2019																																																																																											
Norme armonizzate	Normas harmonizadas	EN IEC 63000:2018	EN IEC 61000-6-3:2021																																																																																											
hawo GmbH Obere Au 2-4 74847 Obrigheim / Germany	T + 49 (0) 6261 / 9770-0 F + 49 (0) 6261 / 9770-69 info@hawo.com www.hawo.com	Amtsgericht Mannheim: HRB 441011 Geschäftsführer: Christian Wolf Firmensitz: Obrigheim	This document and the contents thereof are considered proprietary and confidential information of hawo and disclosure to unauthorized individuals or dissemination, publication, or copying is prohibited without prior written consent by hawo GmbH, 74847 Obrigheim, Germany.																																																																																											

9 910 007, Member 2.01

hawo

hawo GmbH
Obere Au 2 – 4
74847 Obrigheim
Germany
T +49 (0) 6261 / 9770 - 0
F +49 (0) 6261 / 9770 - 69
info@hawo.com
www.hawo.com

hawo USA
321 North Clark Street
Suite 1425
Chicago IL 60654
T +1 312 585 8329
F +1 312 644 0738
info@hawo-usa.com
www.hawo-usa.com

hawo ASIA
25 International Business Park
#03-105 German Centre
Singapore 609916
T +65 6433 5339
F +65 6433 5359
info@hawo-asia.com
www.hawo-asia.com