
White Edition

maXium®

Innovative technology in white!

Where every look and every hand count,
straightforward and easy operation is top priority.
The maXium® unit offers you both.

Reduction is perfection Highest level of simplicity

Nowadays, a trend towards “simplified dialog”, i.e. improved communications with the user, is clearly underway in all fields of medical engineering. Especially during extensive surgeries, extremely complex background processes are going on. Particularly in this field, easy and at the same time intuitive handling of the devices is therefore increasingly gaining significance. The surgeon and his team want to have at their fingertips just those functions and information that they really need for their current tasks, especially in stressful situations.

The maXium® meets these high demands. With this future-oriented generation of HF units, KLS Martin sets entirely new standards in terms of ease of use and performance. maXium® – this is the innovative symbiosis of state-of-the-art high-frequency technology and a whole range of user-friendly assets. This is exactly what makes this unit a “high-end reference system” in the field of multi-functional HF surgery.


maXium®

The unit for all applications

Nowadays, the technical requirements in electrosurgery are particularly complex and diverse. Innovative technology as well as safe and easy handling combined with utmost efficiency are among the most important requirements that have to be met by surgical units.

This is exactly what maXium® offers: The unit exhibits a wide range of precisely regulated currents, enabling you to choose between a great number of applications, such a monopolar and bipolar cutting and coagulation. Besides the standard current types PureCut and BlendCut, you have additional, application-specific currents at your disposal.

But it has still more to offer: The maXium® provides you with important safety mechanisms such as spark control and characteristics control, benefits that guarantee you optimum power adjustment at any time.


Spark control

- Automatic power adjustment to the necessary minimum and thus reduced spark formation – patient-friendly application.
- Constant cutting quality, irrespective of the tissue properties and the electrode geometry

NE indicator


Using the PCS system (Patient Control System – a monitoring system for split neutral electrodes) and the NE indicator, not only the status of the applied neutral electrode but also the application quality is displayed - continuously throughout the entire application.

NE baby

When using small-area neutral electrodes, undesired temperature rises underneath the neutral electrode may occur, depending on the current used. Using the NE baby function, unsuitable current types are blocked and the maximum output power of permissible current types is limited.

Intuitive. Fast. Straightforward. The **maXium®** and its operability


Control takes just a turn

Thanks to its ingenious “Quick Step” control, the maXium® can be easily operated with one hand in a split second. Due to its multifunctional rotary switch, all parameters can be set and changed any time according to your individual needs. The patented backlight additionally indicates the respective status of the rotary switch.

The system’s memory offers you a multitude of application possibilities through 99 storage positions for user-defined programs, plus a preset standard program. These can be stored under the respective medical specialization, the specific indication or a freely definable name – just as you prefer. The list of favorites enables you to access the “last used programs” with ease and call them up in a second when you need them.

Everything at a glance: The maxi-display

Whether it's monopolar, bipolar, cutting or coagulation: The type of the active program, the set output power, the selected type of current and the chosen switching mode (hand or foot switch) can be checked at a glance in the maxi-display. In doing so, the various settings are displayed simultaneously for all four working channels. The unit does not only show you the type of the neutral electrode used via an icon - the NE indicator also continuously provides you information about the current application status on the patient.

Available in three versions

As a standard feature, the maXium® offers you four working channels, each of them with two monopolar and two bipolar outputs. These are available in three different configurations:

- **m-version:** with socket module for KLS Martin accessories
- **i-version:** with socket module for international accessories
- **e-version:** with socket module for Erbe accessories

These three different socket modules allow for the connection of plug-compatible third-party accessories to the maXium® unit without any problems.

Expandable through three port expanders


Each of the four outputs can be selected individually, either via the channel selector button or the rotary switch. The bipolar outputs can additionally be expanded by using bipolar port expanders. This enables you to connect and use a total of four bipolar instruments – without using a second HF unit.

The combination with the bipolar expanders makes it possible to connect e.g. bipolar scissors, bipolar forceps and bipolar sealing instruments simultaneously and configure them independently from each other during extensive surgical interventions. For combining the instruments with KLS Martin, international or Erbe plugs, you have the respective appropriate expanders at your disposal.

Port expanders for instruments:

– with KLS Martin plug

– with international plug

– with Erbe plug

As flexible as you are
The **maXium**® fits in any environment

core nova

Integrated OR systems

Integrated surgical systems are increasingly gaining significance. KLS Martin takes account of this future trend: Now the maxium® unit can be combined with the system CORE nova® of Richard Wolf Company. This makes it possible to change all settings and use all programs required individually for the application out of the sterile area.

Solutions for hand and foot – for those who are quick on the uptake

Program switching with the SWAP mode

The fastest way of changing settings is the SWAP mode. Using this function, you are provided the opportunity to change between two preset programs via the foot switch or the push-buttons on the electrode handle.

The advantages of the SWAP mode:

- Switching directly by the surgeon (no non-sterile person is needed)
- Use of different settings during the surgery
- Optimum performance for all work steps
- Easy handling for creation of SWAP programs
- Switching via:
 - the convenient foot switch with additional button
 - any existing handle with two push-buttons

Plug & Play – automatic identification of instruments with intelligent IQ

With hardware version 06 the name of the connected instrument is automatically displayed by maxium® as soon as an IQ instrument is plugged in, and the current associated with the instrument is set to a matching default value – it simply could not be more straightforward!


If you do not wish to work with the instrument-specific default values, simply change them as required and save the new values in the usual, program-specific manner. The advantage of this method is that when the IQ instrument is plugged in again, your customized values will appear instead of the default values.

The following items are currently available with IQ functionality together with the SealSafe® IQ current mode:

From software version V3.400: marClamp® IQ, marSeal IQ, marSeal Slim IQ and marSeal 5 IQ

From software version V3.412: marClamp® Cut IQ

From software version V3.423: marSeal 5 plus


Displaying the name of the connected instrument

Use of the two bipolar connecting sockets

Automatic assignment of the required current type and set-up of the required performance parameters

Bipolar resection


✓ Effective ✓ Compatible ✓ Safe

With optimized current types and booster pulse

In contrast to classic, monopolar resection, where an electrically non-conductive irrigation liquid is used, the bipolar resection technique exclusively uses an electrically conductive saline solution (NaCl) for irrigation. This provides two advantages for the surgeon: First, the risk of a TUR syndrome is reduced; second, the duration of the surgery can be extended.

Special types of current have been developed specifically for bipolar resection. They feature a booster pulse for starting the cutting process, making it possible to use large cutting loops. Since they exhibit a high resection rate, they are particularly advantageous for TUR-P applications. Using the booster pulse, the cutting process can be initiated quickly and without compressing the tissue with the wire loop.

Examples of bipolar resection

The bipolar resection option gives you the opportunity to use the maXium® unit for bipolar cutting and coagulation in transurethral/transcervical interventions as well. This operating mode requires the use of a bipolar resectoscope manufactured by Storz, Wolf or Olympus.

Main application areas


- Transurethral resection (TUR-B, TUR-P)
- Transcervical resection (TCR)

General advantages of bipolar resection

- Resection can be performed with saline solution
- No risk of causing a TUR syndrome
- No time limits for the intervention
- No neutral electrode required

Special advantages of bipolar resection by KLS Martin


- Excellent cutting effect due to “booster” impulse
- Larger loops possible for faster resection
- Longer service life of the loop electrodes
- No warmed saline solution necessary

Intelligent vessel sealing with
SealSafe® IQ, **marSeal IQ**,
marClamp® Cut IQ and **marClamp® IQ**

Up to 25 % faster sealing!

By optimizing the sealing process, a 25% shorter sealing time is achieved.

Vessel sealing up to ø 7 mm
Lateral thermal expansion < 2 mm
Burst pressure at least 360 mmHg

Safe – vessel sealing with SealSafe® IQ

Bipolar vessel sealing systems allow for sealing vessels and tissue bundles effectively and permanently without any previous dissection or detailed exposure. With the maXium® HF generator, the SealSafe® IQ current, the marClamp® IQ and the marSeal IQ, KLS Martin offers you a fast, reliable and economical vessel sealing system.

Thanks to the precise adaptation of the SealSafe® IQ current to the applications and instruments, exclusively the tissue located between the instrument's jaws is sealed – lateral thermal tissue damage and tissue adherence can thus be kept to an absolute minimum. Besides, the SealSafe® IQ current has been optimized to prevent tissue carbonization to full extent. Due to the effective combination of the constant high pressure applied to the tissue and the high electric current, a permanently sealed zone is reliably created on arteries, veins or tissue bundles.

marSeal Slim IQ

marSeal IQ

marSeal 5 IQ

marSeal 5 plus

marSeal 5 plus

marSeal IQ – one instrument for two work steps

Combined with the maXium® unit and the optimally adapted SealSafe® IQ current, vessels and tissue bundles can be effectively sealed and tissue damage caused by lateral heat propagation can largely be avoided. Adhesion of tissue to the contact surfaces can also be reduced to a minimum. The required constant pressure on the tissue to be sealed is ensured by a ratchet integrated into the handle.

The marSeal IQ instruments provide the optimum basis for easy and tissue-saving sealing with the possibility of subsequent dissection – without any need for exchanging instruments. Thanks to its modular design offering different shaft lengths, the marSeal IQ can be used for laparoscopic applications as well as open surgery. The additional rotation of the tubular shaft by 360 degrees enables easy application of the instrument and fatigue-free working. Since the marSeal can be completely taken apart, it is additionally easy to clean and autoclavable at 134°C (273°F).


The following marSeal IQ instruments are currently available: marSeal IQ, marSeal Slim IQ, marSeal 5 IQ and marSeal 5 plus. For further details, please order our separate marSeal brochure marSeal and marSeal 5 plus!

Application of marSeal IQ in gynaecology:

- Laparoscopic Supracervical Hysterectomy (LASH)
- Laparoscopically Assisted Vaginal Hysterectomy (LAVH)
- Total Laparoscopic Hysterectomy (TLH)
- Vaginal hysterectomy
- Abdominal hysterectomy
- Ovariectomy and salpingectomy
- Adhesiolysis
- Cystectomy

Application of marSeal IQ in general surgery:

- Colectomy
- Nissen Fundoplicatio
- Gastrectomy
- Splenectomy
- Adrenalectomy
- Nephrectomy
- Cholecystectomy

marClamp® Cut IQ

marClamp® Cut IQ instruments and maXium® – the sealing system for fast and safe working

The reusable sealing instrument marClamp® Cut IQ enables tissue preparation and the permanent sealing of veins, arteries and tissue layers with subsequent separation – **without any change of instrument!**

The switch integrated in the instrument automatically triggers the special SealSafe® IQ sealing current once the tissue has been gripped, so no additional footswitch is required.

When sealing has been performed, the current is stopped and an acoustic signal sounds.

The sterile disposable blade guarantees an optimal cutting result at all times.

For more information, please request our separate marClamp® Cut IQ brochure!

Applications of marClamp® Cut IQ

ENT

- Thyroidectomy
- Neck preparation
- Laryngectomy

Gynaecology

- Lymphadenectomy
- Adhesiolysis
- Cystectomy

Visceral surgery

- Colectomy
- Gastrectomy
- Splenectomy
- Adrenalectomy
- Nephrectomy

Urology

- Cystectomy
- Hydrocele
- Partial cystectomy

The appropriate accessories –
the basis for a successful surgery

Tissue-friendly and precise – the marCut® bipolar scissors

Even more safety for patients: Combined with the maXium® unit and the Forfex current, the marCut® bipolar scissors provide you with a system that enables the precise, bloodless and rapid dissection of tissue. Since only the tissue that is located between the scissor blades is integrated into the electric circuit, the energy consumed during marCut® applications is significantly reduced compared to the monopolar procedure. With its special surface hardening and its unilateral ceramic coating, the durable scissor blade provides for a consistently sharp cutting edge, which can be re-sharpened several times. The marCut® is available in standard, slim line and super slim line versions – each of them in different lengths as well.

nonStick red – bipolar forceps with non-stick effect

Combined with the MicroCoag and MacroCoag currents, the new nonStick red enables you to perform bipolar coagulation in all open-surgical applications – without tissue adherence, permanently and reliably. There is a multitude of various models available for this purpose.

Please order our separate brochures marCut® and nonStick red!

maXium®, maXium® Beamer and maXium® smart Cart. Together they form the first-class system for open and endoscopic argon plasma surgery.

maXium®

maXium® Beamer

When performing surgeries with the maXium® Beamer, tissue carbonization is reduced and the wound healing process thus accelerated. Due to the low penetration depth, the risk of perforations is significantly reduced as well.


maXium® smart Vac

maXium® smart Cart

The maXium® smart Cart exhibits a compact and ergonomic design to securely accommodate the entire equipment.

We provide the technology –
to put you in control

ARGON
↓
6.0 l/min

ARGON
Endo
0.5 l/min

ARGON
⏏
1.0 l/min

ARGON
Beam
3.0 l/min

Fast and gentle: Argon makes it possible

In argon plasma surgery, argon gas is ionized with high-frequency current. This produces conductive plasma. The result is extremely homogeneous, flexible and superficial coagulation. During application, the argon gas is always directed at tissue portions that are still bleeding and not yet coagulated.

This process facilitates rapid hemostasis, which is of particular importance especially when dealing with diffuse hemorrhages.

Owing to this, you can significantly reduce blood loss and shorten the operating time. Fractionated currents with low voltage are available specifically for endoscopic applications. They offer you the highest degree of safety and efficiency in interventions with gastrointestinal and tracheobronchial probes.

For further details, please order our separate MABS brochure!

The perfect surgical team:
maXium® and **marVac®**

smartLap

Electrodes with smoke evacuation
for laparoscopy

Various clip-on applicators
combined with the marVac®
smoke evacuator provide reliable
protection from smoke gases.

**Reliable protection from particles and viruses in the operating room:
the marVac® smoke evacuator**

Thanks to its special filter technology, the marVac® eliminates viruses and particles particularly effectively. The suction process is initiated simultaneously with the activation of the HF surgical unit.

Individual suction parameters can be called up and saved directly in the user program of the maXium® HF unit. Thanks to the different connection diameters available, various tubes can be connected without using adapters.

Furthermore, the marVac® impresses with its TURBO and LAP modes as well as the possibility to connect it to a central OR suction system. The maXium® unit and the marVac® smoke evacuator can be accommodated on the maXium® Cart or alternatively integrated into ceiling mounts.

For further details, please order our marVac® brochure!

TURBO mode

Huge quantities of smoke can be eliminated effectively from the operating room.

LAP mode

Automatic reduction of maximum suction performance during laparoscopic interventions

Active suction/ Standby suction

Active suction, standby suction and delay time can be selected as required before and during the surgical procedure. These parameters can be easily adjusted by pressing the up and down buttons.

Modification of parameters

Six good reasons for the **maXium®**

Easy handling – perfect result

The “Quick Step” control integrated into the maXium® unit allows for particularly easy handling - at all times and for all types of applications.

Plug & Play with IQ

The maXium® unit provides you the opportunity to connect encoded KLS Martin IQ instruments. There is a multitude of special instruments available for this purpose.

Maxi-display – maximum sight

The extra-large display immediately and clearly shows you all active parameters.

Always the appropriate current type at hand

There are numerous different current types available for all applications, among them fractionated currents and currents for bipolar vessel sealing.

Perfect cutting results

Using the dynamic characteristics control and the accurate spark control, you obtain cutting currents of up to 360 watts and optimum cutting results.

Preselection of medical discipline

Using the medical discipline preselection function, the maXium® is reduced to the current types and programs designed for the specific medical discipline. Focus on what's essential.

Technical data maXium®

Mains connection	
Power requirements	100 V – 240 V ± 10%
Mains current	max. 6.3 A
Mains fuse	6.3 A
Nominal frequencies	300/400/600 kHz
Power input in standby mode	40 Watt
Max. power input	600 Watt

Output power

Cutting power	max. 360 Watt
Coagulation power	max. 320 Watt

Safety

PCS and NE indicator	
Error log	
Activation log	

Weight and dimensions

Weight	8.3 kg
Width x height x depth	390 x 182 x 436 mm

Standards

Classification	II b
Protection class	I acc. to DIN EN 60601-1
Type of applied part	CF; defibrillation-proof
Approval	CE 0297

Technical data maXium® smart Cart

4 equipotential bonding connectors	
4 dual castors with 2 brakes	

Weight and dimensions

Weight	40 kg
Width x height x depth	755 x 1080 x 655 mm

The requirements for integrating into CORE nova®

• CORE nova®	
• Required hardware by Richard Wolf:	
- Hardware CORE.CONNECT	
- Software driver maXium®	
- Serial data cable 1 m	
• maXium® with hardware 06 and software from V3.398	
• maXium® option CORE nova®	

Technical data maXium® Beamer

Power supply	via maXium® unit
Device fuse	M 1.6 A
Power input	max. 40 Watt
Argon gas flow CUT	0.1 to 12 l/min ± 20% 0.0 l/min = off
Argon gas flow COAG	0.1 to 12 l/min ± 20 % 0.0 l/min = off
LF leakage currents	complying with IEC 601, Part 1 (testing in connection with maxium® HF generator)
Type of applied part	CF; defibrillation-proof
Duty type	INT 10 s/30 s (= duty factor of 25%)

Weight and dimensions

Weight	3.5 kg
Width	390 mm
Height (w/o ball-head bolts)	108 mm
Depth	425 mm

Standards

Classification	II b
Protection class	I acc. to DIN EN 60601-1
Type of applied part	CF; defibrillation-proof
Approval	CE 0297

Description

Description	Item No.
maXium®, m-version	80-042-00-04
maXium®, i-version	80-042-02-04
maXium®, e-version	80-042-04-04
maXium® Beamer	80-044-00-04
maXium® TUR	80-093-00-04
maXium® CORE nova®	80-093-02-04
maXium® smart VAC, with LAP mode, 220 V	80-062-00-04
maXium® smart VAC, 220 V	80-062-01-04
MCB connecting cable	80-091-01-04
maXium® smart Cart	80-048-00-04
Standard rail for maXium® smart cart	80-048-01-04
Separating transformer for maXium® smart cart	80-048-02-04
Device holder for maXium® white edition	80-048-04-04
Gas bottle holder for maXium® smart cart	80-048-05-04
Mounting kit for maXium® smart Vac	80-048-06-04

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co. Ltd.

Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Nippon Martin K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 505 7838
malaysia@klsmartin.com

KLS Martin de México S.A. de C.V.

Mexico City · Mexico
Tel. +52 55 7572 0944
mexico@klsmartin.com

Martin Nederland/Marned B.V.

Huizen · The Netherlands
Tel. +31 35 523 45 38
nederland@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

KLS Martin Taiwan Ltd.

Taipei 106 · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

Reading · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

Do you know how to get all important information
about the KLS Martin Energy Devices?

Please download KLS Martin App Energy Devices!
The App is available for Android and iOS.

[https://itunes.apple.com/de/app/
klm-martin-energy-devices/id1198171415?l=de&ls=1&mt=8](https://itunes.apple.com/de/app/klm-martin-energy-devices/id1198171415?l=de&ls=1&mt=8)

[https://play.google.com/store/apps/
details?id=com.klsmartin.energydevices](https://play.google.com/store/apps/details?id=com.klsmartin.energydevices)

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
P.O. Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com


Rigid and flexible probes
for open and endoscopic argon surgery

Rigid and flexible Probes for open and endoscopic Argon Surgery

Argon-supported HF surgery has been offering advantages for surface coagulation for many years. The application fields of argon-supported HF surgery include:


- open surgery
- laparoscopic surgery
- pneumology
- gastroenterology


Working principle

By ionizing the argon gas atmosphere using HF current (high-voltage spray coagulation current), electrically conductive argon plasma is generated.

Electrodes of different shape and length are available to users, thus providing for cutting as well as large-area coagulation.

Advantages of argon surgery:

- Reliable coagulation with minimal traumatization of the tissue or organ
- Low blood loss and shortened operating times
- Minimal vaporization at a very low penetration depth (0.5 to max. 3 mm), therefore significantly reduced perforation risk
- No carbonization, therefore faster wound healing
- Significantly reduced smoke development, improved view of the surgical site
- Non-contact method – no tissue sticking to electrode
- Less bronchial and pleural fistula formation on lung parenchyma

Application examples for Argon Probes

Argon Cut

Argon Beam

Argon Endo

**Pulsed Argon
Slow Rep.**

**Pulsed Argon
Fast Rep.**


Open surgery examples:

- General surgery
 - Partial liver resection, colon resection
- ENT
 - Tonsillectomy
- Cardiac/thoracic surgery
 - Median sternotomy
- Emergency surgery
 - Large-surface hemostasis during skin grafting procedures
- Plastic surgery
 - Mammary reduction

Endoscopic surgery examples:

- Gastroenterology
 - Gastric ulcer
 - Esophagus
 - Angiodysplasias
 - Hemostasis of diffuse bleeding
 - Conditioning prior to fistula sealing
 - Hemostasis after polypectomy
- Gynecology
 - Hemostasis after myomectomy
- Laparoscopy
 - Hemostasis after cholecystectomy
- Pneumology
 - Coagulation of superficial hemorrhages
 - Coagulation in trachea, primary/secondary bronchi

MABS Snare probes — the two-in-one solution

MABS Snare probes are combination probes allowing you to remove polyps with subsequent argon coagulation of the site.


The great advantage of these probes: you don't need to exchange instruments for performing the coagulation in a second step. Following ablation, you just pull the snare back into the instrument to use the probe as a "normal" MABS beam probe. Moreover, the probe can be conveniently used in this way for coagulating/vaporizing minor polyps.

Why coagulation following polypectomy?
The answer is that additional coagulation cuts the recurrence rate.

The ideal combination between resection and argon beamer coagulation

Step 1: resection

Step 2: argon beamer coagulation

Flexible MABS probes for gastroenterological and endobronchial applications

Flexible MABS electrodes have the following features in common:

- distal ceramic nozzle
- autoclavable at 134°C / 273°F (only reusable probes)
- reduced gas consumption (50% less than previous models)

80-181-31-04

MABS rinsing adapter for reusable argon probes,
5 item(s)/Pack

80-181-30-04

MABS connecting cable for flexible probes (disposable and reusable),
cable length 2.5 m, for HF current and argon gas,
HF-current and gas-flow activation via foot switch,
autoclavable at 134°C / 273°F (only reusable probes)


Argon probe (reusable)

Argon probe (disposable)

MABS probes


Item No.	Designation		Diameter	Length	Item(s)/Pack
80-181-23-04	MABS GIT probe, reusable		2.3 mm	2.3 m	1
80-181-24-04	MABS GIT probe, reusable		3.2 mm	2.3 m	1
80-181-25-04	MABS TBS probe, disposable		1.8 mm	1.6 m	10
80-181-26-04	MABS GIT probe, disposable		1.8 mm	3.2 m	10
80-181-27-04	MABS GIT probe, disposable		2.3 mm	2.3 m	10
80-181-28-04	MABS GIT probe, disposable		3.2 mm	2.3 m	10
80-181-29-04	MABS GIT probe, disposable		2.3 mm	3.4 m	10
80-181-32-04	MABS GIT probe, Side Fire, disposable		2.3 mm	2.3 m	10

1/1

80-181-41-04

MABS Snare probe (multifilament, disposable)

80-289-40-04

4 m/12 ft.

Monopolar connecting cable
for KLS Martin HF units

80-289-41-04

4 m/12 ft.

Monopolar connecting cable
for HF units with 3-pin connectors

80-289-42-04

4 m/12 ft.

Monopolar connecting cable for maxium® "e"
version and Erbe ICC/ACC/VIO units

MABS snare probes

Item No.	Designation		Diameter	Length	Items/Pack
80-181-41-04	MABS GIT snare	15 mm, multifil., disposable	2.5 mm	2.3 m	5
80-181-42-04	MABS GIT snare	30 mm, multifil., disposable	2.5 mm	2.3 m	5
80-181-43-04	MABS GIT snare	15 mm, monofil., disposable	2.5 mm	2.3 m	5
80-181-44-04	MABS GIT snare	30 mm, monofil., disposable	2.5 mm	2.3 m	5

monofil

multifil

Rigid and flexible MABS probes for open and endoscopic applications

Rigid MABS electrodes have the following features in common:

- insulated, rigid shaft with a diameter of 5 mm
- distal ceramic nozzle
- autoclavable at 134°C (273°F)

80-181-02-04

MABS handle for rigid applicators
Two pushbuttons for coagulation and cutting
Connecting cable, 4.5 m/15 ft., for HF current and argon gas
Autoclavable at 134°C (273°F)

80-181-07-04

MABS needle electrode, axially adjustable

80-181-08-04

MABS needle electrode, axially adjustable

80-181-05-04

Fixing cap, suitable for all rigid electrodes

80-181-10-04

MABS beam electrode for open surgery,
with high-temperature-proof ignition tip


80-181-11-04

MABS beam electrode for open surgery,
with high-temperature-proof ignition tip

80-181-12-04

MABS beam electrode for laparoscopy and pelviscopy, with high-temperature-proof ignition tip and fixing cap


80-181-13-04

MABS lancet electrode for open surgery

80-181-14-04

MABS lancet electrode for open surgery

80-181-15-04

MABS needle electrode for open surgery


80-181-16-04

MABS needle electrode for open surgery

80-181-60-04

MABS beam electrode, 320 mm, flexible

80-181-62-04

MABS beam electrode, 115 mm, flexible

KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co., Ltd

Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Nippon Martin K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 506 2380
malaysia@klsmartin.com

KLS Martin de México S.A. de C.V.

Mexico City · Mexico
Tel. +52 55 7572 0944
mexico@klsmartin.com

Martin Nederland/Marned B.V.

Huizen · Netherlands
Tel. +31 35 523 45 38
infonl@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Moscow · Russia
Tel. +7 499 792 76 19
russia@klsmartin.com

KLS Martin Taiwan Ltd.

Taipei 106 · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

Reading · United Kingdom
Tel. +44 118 467 1500
uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

Do you know how to get all important information
about the KLS Martin Energy Devices?

Please download KLS Martin App Energy Devices!
The App is available for Android and iOS.

[https://itunes.apple.com/de/app/
klm-martin-energy-devices/id1198171415?l=de&ls=1&mt=8](https://itunes.apple.com/de/app/klm-martin-energy-devices/id1198171415?l=de&ls=1&mt=8)

[https://play.google.com/store/apps/
details?id=com.klsmartin.energydevices](https://play.google.com/store/apps/details?id=com.klsmartin.energydevices)

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Option integration into Karl Storz OR1 NEO system

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Option integration into CORE nova®/Richard Wolf system

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For use only with the maXium® generator

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Double-pedal foot switch, explosion- and anesthetics-proof (AP-tested), IP X8, immersible

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maXium® double-pedal foot switch with SWAP mode Eco, AP-tested (anesthetics-proof), IP X8, immersible

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Foot switch, explosion- and anesthetics-proof (AP-tested), IP X8, immersible

Interruptor de pie, antiexplosivo, controlada su aptitud para anestésicos, IP X8, eintauchbar

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Cable length for all foot switches: 5 m

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Third-party devices with this type
of connector

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

Electrode handle with two pushbuttons, 3-pin plug and permanent cord of 4 m

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Electrode handle with rocker switch, 3-pin plug and permanent cord of 4 m

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Electrode handle with two pushbuttons, with detachable cord

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


Twin-Pad II neutral electrode with two contact surfaces (split), contact surface 200 cm²
(1 x 80-294-40-04/80-294-43-04/80-294-44-04 + 1 x 80-372-00-04 required)

Placa neutrale Twin-Pad II con dos superficies, superficie de contacto de 200 cm²
(preciso 1 x 80-294-40-04/80-294-43-04/80-294-44-04 + 1 x 80-372-00-04)

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When using the Twin Pad II, the use of a compression Class 1 support stocking is recommended.

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Textile band with velcro strap, 90 cm,
washable at 60°C

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


* incl. 80-630-17-04, 80-630-18-04
** incl. 80-630-17-04
Blades have to be ordered separately

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For further information please refer to our separate product brochure.

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* incl. 80-633-12-04, 80-633-13-04

** incl. 80-633-12-04
Blades have to be ordered separately

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Gebrüder Martin GmbH & Co. KG

KLS Martin Platz 1
78532 Tuttlingen
Germany

has implemented and maintains a **Quality Management System**.

Scope:

Development, manufacturing and distribution of medical products, surgical instruments and devices, non-active implants, RF surgery units, operation theatre light systems, ceiling-mounted surgery stands, surgical laser devices and container systems for sterilization and accessories for the above mentioned products.

Through an audit, documented in a report, it was verified that the management system fulfills the requirements of the following standard:

ISO 9001 : 2015

Certificate registration no.	210299 QM15
Certificate unique ID	170738717
Effective date	2020-07-29
Expiry date	2023-07-28
Frankfurt am Main	2020-07-29

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Sigrid Uhlemann
Managing Director

August-Schanz-Straße 21, 60433 Frankfurt am Main,
Tel. +49 (0) 69 95427-300, medical.devices@dqs-med.de

CERTIFICATE

This is to certify that the company

KLS martin
GROUP

Gebrüder Martin GmbH & Co. KG

KLS Martin Platz 1
78532 Tuttlingen
Germany

has implemented and maintains a **Quality Management System**.

Scope:

Development, manufacturing and distribution of medical products, surgical instruments and devices, non-active implants, RF surgery units, operation theatre light supply units, ceiling-mounted supply systems, surgical laser devices and container systems for sterilization and accessories for the above mentioned products.

Through an audit, documented in a report, performed by DQS Medizinprodukte GmbH, it was verified that the management system fulfills the requirements of the following standard:

EN ISO 13485 : 2016

Certificate registration no.	210299 MP2016
Certificate unique ID	170738702
Effective date	2019-08-29
Expiry date	2022-08-28
Frankfurt am Main	2020-03-06

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Sigrid Uhlemann
Managing Director

Dr. Thomas Feldmann
Head of Certification Body

August-Schanz-Straße 21, 60433 Frankfurt am Main,
Tel. +49 (0) 69 95427-300, medical.devices@dqs-med.de

EC-CERTIFICATE

(Full quality assurance system)

This is to certify that the company

KLS martin
GROUP

Gebrüder Martin GmbH & Co. KG

KLS Martin Platz 1
78532 Tuttlingen
Germany

has implemented and maintains a full quality assurance system which applies to the products at every stage from design to final controls.

Through an audit, documented in a report, performed by DQS Medizinprodukte GmbH, it was verified that the management system fulfills the requirements of

Annex II – excluding Section 4 of Council Directive 93/42/EEC concerning medical devices

with respect to the following medical devices:

Active and non active medical devices, surgical instruments and units, implants, medical devices for RF- and laser surgery, Smoke evacuation systems, Accessories for Smoke evacuation systems, ceiling-mounted supply systems and accessories according to annex.

The manufacturer is subject to surveillance according to Annex II, Section 5. The CE marking with the Notified Body Identification Number (0297) may be affixed on the devices listed in the certificate. An EC Design Examination Certificate according to Annex II, Section 4 is required for class III devices covered by this certificate. The certificate is in the case of class I(s) devices (I(s) = class I products placed on the market in sterile conditions) limited to the aspects of manufacture concerned with securing and maintaining sterile conditions. The certificate is in the case of class I(m) devices (I(m) = class I devices with a measuring function) limited to the aspects of manufacture concerned with the conformity of the products with the metrological requirements.

Certificate registration no.	210299 MR2
Certificate unique ID	170763177
Effective date	2020-01-13
Expiry date	2024-05-26
Frankfurt am Main	2020-01-13

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Sigrid Uhlemann
Managing Director

Dr. Thomas Feldmann
Head of Certification Body

August-Schanz-Straße 21, 60433 Frankfurt am Main,
Tel. +49 (0) 69 95427-300, medical.devices@dqs-med.de

DQS Medizinprodukte GmbH is a Notified Body according to Council Directive 93/42/EEC concerning medical devices with the Identification Number 0297.

Annex to certificate

Certificate registration No.: 210299 MR2

Certificate unique ID: 170763177

Effective date: 2020-01-13

Gebrüder Martin GmbH & Co. KG

KLS Martin Platz 1
78532 Tuttlingen
Germany

Device family	Devices	Class
RF-surgery units	Medical RF-surgery units for general and special applications	IIb
Accessories for RF-surgery units	Argon Beamer Handles for RF-surgery Electrodes Forceps and scissors Clamps, forceps, NervScout	IIb IIa / IIb
Laser surgery units	surgical laser units for common and particular use	IIb
Accessories for laser surgery units	Laser contact probe Adapter for laser fibres Connection fibres Fibre optics Bipolar suction and irrigation units Micromanipulator	IIa / IIb
Smoke evacuation systems	products for smoke evacuation	IIa
Accessories for Smoke evacuation systems	adapters	IIa
Implants and Instruments	Distractor Trocator tip for activation arm Condylar implant Washers Meshes (for cranioplasty) Bone plates Bone screws Bone nails Fixation screws Fixation screws (lug) UHP (prosthesis for wrist joint) Tantal wire (orthodontic) Bone clips Bone wire Finger joint CranioXpand marPOR Polyethylene Implants	IIa / IIb

Annex to certificate

Certificate registration No.: 210299 MR2

Certificate unique ID: 170763177

Effective date: 2020-01-13

Gebrüder Martin GmbH & Co. KG

KLS Martin Platz 1
78532 Tuttlingen
Germany

Device Family	Device	Class
Resorbable Implants	Resorb X meshes	III
	Resorb X plates	
	Resorb X screws	
	Resorb X Pins	
Surgical instruments and units	Sterile scalpels	Ila
	Sterile blades	Ila
	Guide pins	Ilb
	Guide wires	Ilb
	Drill bits OMF surgery	Ila
	BOS driver	Ila
	SonicWeld Rx (basis device, handpieces)	Ila
	SonicWeld Rx (sonotrodes)	Ila
	Drill bits / Saw blades orthopaedic	Ila
	BOS-Battery pack sterile	Is
	Bending Templates	Is
Pendants	Service Head	Ilb
Wound spreading systems	marTract	Ila
Surgical Instruments and units	Suction Units, Suction Cannulae, Irrigation Cannulae	Ila