

INgezim BLV Confirmation

R.12.BLV.K1

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

INgezim BLV Confirmation

R.12.BLV.K1

INgezim BLV Confirmation is an immunoenzymatic assay based on an indirect ELISA technique, which uses a monoclonal antibody (MAb) specific of bovine immunoglobulins. This assay uses plates coated with positive antigen (inactivated BLV) and negative antigen.

TECHNICAL BASIS OF THE KIT

1. Plates are coated with inactivated BLV antigen and negative antigen respectively. Samples are added to both plates and incubated.
2. If the sample contains specific antibodies to BLV, they will bind to the antigen.
3. When a MAb-PO specific for bovine immunoglobulins is added, it will bind to the Igs of the sample previously bound to the antigen. This binding is detected by the development of a colorimetric reaction after the substrate addition.

APPLICATION

Confirmatory assay for detection of specific antibodies to Bovine Leukaemia Virus (BLV) in bovine serum samples.

INTERPRETATION OF THE RESULTS

One cut off is used for results interpretation. Samples will be considered **Positive** if their OD value corrected with respect to the OD value obtained with negative antigen is higher than the cut off. Samples will be considered **Negative** if this corrected OD value is lower than the cut off.

VALIDATION

1. Use of field sera.

A study was made using a set of 400 sera previously classified by INGEZIM BLV COMPAC:

- 29 sera positive from Canada and Europe
- 371 sera negative from Holland, Spain, Poland and Canada

Classified	n	Pos	Doub	Neg
Negative	371	4	0	367
Positive	29	28	0	1

Correspondence between both assays was 98,7%

2. Use of O.I.E. reference sera

In order to verify the sensitivity of the assay, E-4 and E05 Community Reference Sera were requested, and, following the O.I.E. legislation, they were assayed at a 1/10 dilution in negative serum (as weak positive serum) and diluted 1/100 in negative sera. Positive results were obtained in all cases. The OIE Negative Reference Serum showed negative results in the same assay. Therefore, we conclude that INGEZIM BLV Confirmation maintains the level of sensitivity and specificity required by the O.I.E legislation, and therefore can be used as a confirmation assay.

COMPOSITION OF THE KIT

- Microtitration plates coated with BLV antigen.
- Microtitration plates coated with negative antigen.
- Vial with Positive Control.
- Vial with Negative Control.
- Vial with Conjugate.
- Bottles with Washing Solution concentrated.
- Bottles with diluent
- Bottle with substrate (TMB) ready to use.
- Bottle with stop solution.

PRODUCT MANUFACTURED BY INGENASA
REGISTRATION No: 3471 RD

Shelf life: 18 months

Ed.020217

INGENASA

Store at 2°C-8°C

Ed.020217