## Copyrights and Disclaimers

©2016, Immucor, Inc. NEO Iris™ is a trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrol-de@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged.

No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied

by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

### **Document Revision History**

Date	Version	Chapters	Description
MAY 2016	NEO Iris_EU-001-100	N/A	First Version of NEO Iris Operating Manual
JUL 2016	NEO Iris_EU-001-101	Chapter 12: Limitations of Use and Warnings Attachment 1	Update of the limitations on interfering substances for Lipemia and Icterus (Levels defined)
	June 2017	Copyrights and Disclaimers	Manufacturer new address and Document Revision History
	NEO Iris_EU-001-101	Chapter 1: Introduction to the NEO Iris	About the NEO Iris – Intended use Change Instrument label including the CE logo
JUN 2017	NEO Iris_EU-001-102 (A-I)	Attachment 1: NEO Iris Operator Manual	Manufacturer new address and Document Revision History List of assays updated Change the Anti-A positive cutoff range for samples from 70 to 100 to 58 to 100 for ABO assays
	NEO Iris_EU-001-101 (A-II)	Attachment 2: NEO Iris Operator Manual	Manufacturer new address Remove FDA clearance claim
	November 2018	Copyrights and	Update of DRH

Date	Version	Chapters	Description
		Disclaimers	
	NEO Iris_EU-001-101	02-03- TOC_about_ NEO_Iris_EU	Warning included (Hot surface)
NOV 2018	NEO Iris_EU-001-102	Chapter 1: Introduction to the NEO Iris	References included for: - CMOS cameras - Change from CSA to Nemko certification in label
	NEO Iris_EU-001-101	Chapter 2: Harware Components	Reference included: Complementary Metal-Oxide Semiconductor (CMOS)
	NEO Iris_EU-001-101	Chapter 3: System Software Navigation	Autologoff notes included in Logging in section
NEO Iris NEO Iris NEO Iris	NEO Iris_EU-001-101	Chapter 5: Instrument Startup	Autologoff section included
	NEO Iris_EU-001-101	Chapter 10: Maintaining the NEO Iris	Limitation included in Maintenance Overview section
	NEO Iris_EU-001-102	Chapter 12: Limitations of Use and Warnings	Limitation included in: - Instrument Start-up - Maintaining the NEO Iris - Preparing the NEO Iris for the First Use
	NEO Iris_EU-001-101	Appendix A: Preparing the NEO Iris for the First Use	Update in: - Warning and Limitations - Standards
	NEO Iris_EU-001-101	Appendix C: Hardware	Update in General Data and Personal Computer tables

Date	Version	Chapters	Description
		Technical Data	
	NEO Iris_EU-001-101	Glossary of Terms	CCD camera reference removed
	NEO Iris_EU-001-103	Attachment 1	Capture-CMV included
DEC 2020	NEO Iris_EU-001-102	Chapter 10: Maintaining the NEO Iris	Updated recommended cleaning solution in the following maintenance tasks: - Daily: 'Cleaning the Instrument' - Weekly: 'Cleaning the Pipettor Wash Towers' - Weekly: 'Cleaning the Common Waste Container' - Monthly: 'Decontaminate Tubings' including specifications for sodium hypochlorite and steps to prepare solution. Revised instructions for performing the monthly 'Decontaminate Tubings' maintenance task.
MAY 2021	NEO Iris_EU-001-102	Appendix C: Hardware Technical Data	Update of Operating System (to Windows 10) in Personal Computer table
MAR 2023	NEO Iris_EU-001-102	Chapter 12	Update to include reference to reagent Instructions
	NEO Iris_EU-004-103	Attachment 1	for Use related to the anticoagulants

# **Table of Contents**

About This Manual	xi
How This Manual Is Organized	xii
Chapter 1: Introduction to the NEO Iris	1-1
About the NEO Iris	1-2
About the NEO Iris	1-2
Consistent Color Code	1-3
Barcodes	1-6
Reagents	1-10
Racks	1-11
Dynamic Scheduler	1-13
Miscellaneous Information	1-14
Chapter 2: Hardware Components	2-1
The NEO Iris	2-2
Cabinet	2-5
Hood	2-9
Loading Tower	2-11
Transport System	2-14
Plate Carrier	2-17
14-lane and 5-lane Bays	2-18
Pipetting System	2-20
Incubator	2-26
Washer	2-28
Centrifuge	2-31
Camera Reader	2-36
Disposal Information	2-38
Chapter 3: System Software Navigation	
Navigation	
Main Menu Bar	
Machine Monitor	3-59
Status Bar	
Chapter 4: Security	4-1
Assigning Passwords and User Access Rights	
Adding a User	
Editing a User	4-10

Deleting a User4	-12
Default Access Rights4	-14
Changing a Password4	-16
Archive Configuration4	-18
Chapter 5: Instrument Start-Up	5-1
Starting Up	5-2
Logging In and Initialization	5-3
Auto Logoff	5-5
Chapter 6: Instrument Testing Operation	6-1
Using the Start Run Assistant	6-2
Loading Samples	6-3
Downloading Requests from LIS6	-14
Completing the Sample Loading Process6	-16
Loading Reagents and Controls6	-21
Loading Plates6	-27
Starting Processing6	-34
Continuous Loading During Operation6	-36
Chapter 7: Test Results	7-1
Accessing the Results Screen	7-2
Sample View and Plate View Icons and Symbols	7-7
Using Tool Tips	7-9
Viewing Test Details7	-10
Approving Test Results7	-17
Exporting Test Results7	-18
Viewing Archives7	-20
Chapter 8: NEO Iris Reports	8-1
Reports Overview	8-2
Parts of the Report	8-4
Plate Based Reports8	-10
Sample Based Reports8	-13
Current Reports8	-15
Quality Control Reports8	-20
Reagent Reports8	-23
Accessing Plate Based Reports8	-25
Accessing Sample Based Reports8	-27
Printing Reports8	-29
Printing Reagents Report from Test Details8	-35

Test Results and Interpretation	8-37
Chapter 9: System Shutdown	9-1
Logging Out	9-2
Shutting Down the NEO Iris After Operation	
Extended Shutdown of the NEO Iris	9-9
Chapter 10: Maintaining the NEO Iris	
Maintenance Overview	
Daily Maintenance Tasks	
Weekly Maintenance Tasks	
Monthly Maintenance Tasks	
As Needed Maintenance Tasks	
Chapter 11: Troubleshooting the NEO Iris	11-1
The Troubleshooting Process Steps	
Using Error Codes to Troubleshoot	
Troubleshooting Software Failure	11-10
Pipettor Self Check Failures	11-14
Clot Detection Recovery Process	11-18
Troubleshooting Plate Transport Errors	11-19
Troubleshooting Pipettor Errors	11-27
Troubleshooting Centrifuge Errors	11-32
Troubleshooting Incubator Errors	11-44
Troubleshooting Washer Errors	11-50
Troubleshooting Camera Reader Errors	11-58
Troubleshooting 14-lane and 5-lane Bay Errors	11-60
Troubleshooting Plate Tower Errors	11-62
Chapter 12: Limitations of Use and Warnings	
Limitations of Use	
Warnings	12-16
Appendix A: Preparing the NEO Iris for First Use	A-1
Verifying all Parts Are Present	A-2
Environmental Conditions and General Safety Features	A-3
User Safety	A-7
Making the Connections	A-9
Software Installation	A-12
Setting Up the Instrument	A-13
Completing the Post-Installation Check	A-14
Verifying the Installation	A-15

Removal of the InstrumentA-16
Appendix B: Maintenance RecordsB-1
Maintenance FormsB-2
Appendix C: Hardware Technical DataC-1
Hardware Technical DataC-2
Glossary of Terms1
Glossary2
Attachment I: NEO Iris Operator ManualI-1
Copyrights and DisclaimersI-2
Sample RequirementsI-5
Assay DescriptionsI-8
Assay CutoffsI-24
Assay Reagent Component GridI-48
Assay Procedural StepsI-51
Test Results and InterpretationI-82
Attachment II: NEO Iris Operator ManualII-1
Copyrights and DisclaimersII-2
ScopeII-4
Essential Information for CommunicationII-5
Result Message StructureII-10
Host Query Message StructureII-24
Order Message StructureII-27
Message ExamplesII-30
Attachment III: NEO Iris Operator Manual1
Copyrights and Disclaimers2
How this Attachment is Organized4
ABO Titration Assay Description6
ABO Titration Cutoffs and Reagent Components12
ABO Titration Procedural Steps13
ABO Titration Results and Interpretations16
Interface Specification Information20
Attachment IV: NEO Iris Operator Manual1
Copyrights and Disclaimers2
How this Attachment is Organized4
Monoclonal Jka/Jkb Assay Description6
Monoclonal Jka/Jkb Assay Cutoffs and Reagent Components9
Monoclonal Jka/Jkb Assay Procedural Steps11

Monoclonal Jka/Jkb Assay Results and Interpretations	12
Attachment V: NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Modified Antigens Assays Description	6
Modified Antigens Assays Cutoffs and Reagent Components	8
Modified Antigens Assays Procedural Steps	9
Modified Antigens Assays Results and Interpretations	11
Attachment VI: NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Assay Description	6
Cutoffs and Reagent Components	12
Procedural Steps	54
Results and Interpretations	56
Attachment VII: NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Assay Description	6
Cutoffs and Reagent Components	7
Procedural Steps	8
Results and Interpretations	9
Attachment VIII: NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Assay Description	6
Cutoffs and Reagent Components	10
Procedural Steps	32
Results and Interpretations	33
Attachment IX: NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Assay Description	6
Cutoffs and Reagent Components	8
Procedural Steps	9
Results and Interpretations	11
Attachment X: NEO Iris Operator Manual	1

Copyrights and Disclaimers2
How this Attachment is Organized4
Antigen Screening Assay Description5
Antigen Screening Assay Cutoffs7
Antigen Screening Assay Reagent Component Grid7
Antigen Screening Assay Procedural Steps8
Antigen Screening Assay Results and Interpretations
Attachment XI: NEO Iris Operator Manual1
Copyrights and Disclaimers2
How this Attachment is Organized3
Non-ABO Titration Assay Description4
Non-ABO Titration Cutoffs and Reagent Components7
Non-ABO Titration Procedural Steps9
Non-ABO Titration Results and Interpretations
Performance Characteristics12
Interface Specification Information13
Attachment XIV: Capture-R Select based Screening Assays1
Copyrights and Disclaimers2
How this Attachment is Organized4
Description for Assays6
Cutoffs and Reagent Components12
Procedural Steps14
Results and Interpretations16
RedCell Labeling Sheet18
Attachment XV: NEO Iris Operator Manual1
Copyrights and Disclaimers2
How this Attachment is Organized4
Description for ABODFULLH, ABORH_V and Baby_BG3 Assays6
Cutoffs and Reagent Components8
Procedural Steps14
Results and Interpretations16
Attachment XVI: NEO Iris Operator Manual1
Copyrights and Disclaimers2
How this Attachment is Organized4
Description for monoclonal S/s, Fya, Fyb, k Assays6
Monoclonal S/s, Fy <sup>a</sup> , Fy <sup>b</sup> , k Assay Cutoffs and Reagent Components10
Monoclonal S/s, Fya, Fyb, k Assay Procedural Steps

Monoclonal S/s, Fya, Fyb, k Assay Results and Interpretations	17
Interface Specification Information	20
Performance Characteristics	22
Attachment XVIII: NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Description for RHFORMEL and PHENO12	6
Cutoffs and Reagent Components	8
Procedural Steps	14
Results and Interpretations	15
Interface Specification Information	18
Performance Characteristics	20
Attachment XIX: Capture-R Select based Identification Assays	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Description for 8-Cell ID Assays	6
Cutoffs and Reagent Components	11
Procedural Steps	14
Results and Interpretations	17
RedCell Labeling Sheet	19
Attachment XX: NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Description for monoclonal M antigen Assay	6
Monoclonal M Assay Cutoffs and Reagent Components	8
Monoclonal M Assay Procedural Steps	10
Monoclonal M Assay Results and Interpretations	11
Interface Specification Information	14
Performance Characteristics	15
Attachment XXI NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Description for monoclonal N antigen Assay	6
Monoclonal N Assay Cutoffs and Reagent Components	9
Monoclonal N Assay Procedural Steps	11
Monoclonal N Assay Results and Interpretations	12
Interface Specification Information	15

Performance Characteristics16
Attachment XXII: NEO Iris Operator Manual1
Copyrights and Disclaimers2
How this Attachment is Organized4
Description for AG_K (H), PHENO12_2, PHENO16 and PHENO16_2 assays6
Cutoffs and Reagent Components8
Procedural Steps16
Results and Interpretations17
Interface Specification Information20
Performance Characteristics22
Attachment XXIV: NEO Iris Operator Manual1
Copyrights and Disclaimers2
How this Attachment is Organized4
Assay Description6
Cutoffs and Reagent Components11
Procedural Steps44
Results and Interpretations46
Performance Characteristics51
NEO Iris® ASSAY GUIDE:1
Open Channel Anti-D Group Assays1
About this Assay Guide2
Intended Use2
List of Assays and Reagents for the NEO Iris <sup>®</sup> 2
Assay Procedural Steps5
Test Results and Interpretations7
Assay Cutoffs for NEO Iris <sup>®</sup> 9
Reagent Barcodes10
Limitations and Warnings12
Example Labelling Sheet for Open Channel Assays13
Index1
12
52
A2
B2
C3
D4
E5

E CONTRACTOR OF CONT	5
-	
G	5
Н	5
Ι	6
К	6
L	7
М	7
Ν	9
O	9
Ρ	9
Q	11
R	11
S	12
Т	14
U	15
V	15
W	16
Υ	16

# **About This Manual**

#### In This Section

The *NEO Iris Operator Manual* is designed to guide the NEO Iris<sup>™</sup> operator through all procedures required to use and maintain the NEO Iris, including operating procedures, maintenance, and troubleshooting.

This chapter provides high-level information about how this manual is organized.

About This Manual	xi
How This Manual Is Organized	xii

# How This Manual Is Organized

### In This Section

This section describes the organization of this manual, including:

- Format
- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons
- Front Pages
- Chapters

#### Format

This manual is divided into chapters that describe specific aspects of the NEO IRIS structure or functionality.

### **Notational Conventions**

This manual uses a page numbering system that includes a prefix of the chapter number hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current manual version identified using a nine character hyphenated format. The first three characters (NEO) identify the instrument. The second set of three characters identifies this document as the operator manual (001). The final set of three characters identifies the version of the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

If changes are made to individual chapters in between full manual version updates, then chapter version numbers are incremented by single digits. For example, 101 is the first update of a chapter falling between a full manual update from version 1 to version 2.

### Limitations of Use and Warnings

Limitations of use and warnings are located throughout this manual, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings. **Chapter 12 – Limitations of Use and Warnings** contains a complete list of all of the limitations of use and warnings in this manual.

### Use of Icons

Some of the following safety symbols may be displayed on the NEO Iris or may appear in the manual to alert you of warnings or limitations of use, or to direct you to information. Examples are shown below.

Symbol	Type of Warning	Related to
	Laser beam safety warning	Laser beam safety issues
$\sim$	Alternating current	Power supply
	Direct current	Power supply
	Protective conductor terminal	Power supply
	Fuse	Power supply
	On (supply)	Power supply
0	Off (supply)	Power supply
	Limitations and Warnings	Potentially damaging or dangerous outcomes if certain critical procedural steps are ignored or incorrectly executed
4	Warning, risk of electric shock	Potential hazard related to power supply
	Warning, risk of crushing or pinching	Potential hazard resulting in possible injury
	Hot surface warning	Y-Motors at pipettor arms
Ţ	Consult instructions for use	
IVD	In Vitro Diagnostic medical device (IVD)	

	Biological risks	
	Manufacturer	
$\sim$	Manufacturing date	
X	Separate collection for electrical and electronic equipment	
EC REP	Authorized Representative in the European Community	
SN	Serial number	
REF	Catalog number	

## Front Pages

The front pages of this manual include the Copyrights and Disclaimers page and the Table of Contents.

## Chapters

The Table of Contents sequentially lists all chapter contents. The first page of each chapter lists the contents of that chapter.

# **Chapter 1: Introduction to the NEO Iris**

# In This Chapter

This chapter provides an introduction to the NEO Iris.	
CHAPTER 1: INTRODUCTION TO THE NEO IRIS	1-1
About the NEO Iris	1-2
Consistent Color Code	1-3
Barcodes	1-6
Reagents	1-10
Racks	1-11
Dynamic Scheduler	1-13
Miscellaneous Information	1-14

## About the NEO Iris

#### Intended Use

The Neo Iris (NEO Iris) is a microprocessor-controlled instrument to fully automate immunohematology in vitro diagnostic testing of human blood. The NEO Iris automates test processing, result interpretation and data management functions. The NEO Iris is designed to automate standard immunohematology assays using a microplate-based platform. Assays include ABO grouping and Rh (D) typing, detection/identification of IgG red blood cell antibodies, compatibility testing, red blood cell phenotyping, antigen screening and infectious disease screening, such as Cytomegalovirus (CMV).

The NEO Iris is part of the Galileo Family of instruments and is intended for use only with the reagents described in **Attachment 1 for Neo Iris Operator Manual**.



Limitation: The Neo Iris is for in vitro diagnostic use.

### **Principles of Operation**

The NEO Iris is designed to automate standard immunohematology assays and operate as a walk-away system, meaning you can leave the NEO Iris to operate independently for periods of time. This leaves you free to carry out other tasks in the laboratory. Several unified principles have been integrated into the NEO Iris system to support and to simplify the overall system operation.

The NEO Iris is a closed system and can only be used with specified Immucor products.

The NEO Iris is an ergonomically friendly and easy-to-use system. Features of the NEO Iris system have been designed to maximize operator efficiency and thereby minimize result errors.

The NEO Iris is a robotic instrument programmed to move microplates, liquid reagent fluids, and blood sample fluids to different bays and processing areas for a given assay in the correct sequence. Such bays and areas include incubator bays, the microplate washing station, the centrifuge, and the reader.

The NEO Iris plate reader uses CMOS cameras to capture an image of the microplate from underneath. The NEO Iris software calculates a reaction value for each well based on a multi feature image analysis. The NEO Iris then assigns a result and interpretation to the wells based on predefined criteria associated with the calculated reaction value. Some assay protocols require multiple test wells for a given blood sample interpretation, such as ABO and Rh (D) typing.

The NEO Iris uses software to drive its mechanics and data processing. The operator uses hardware in combination with the software to operate and maintain the NEO Iris.

## Consistent Color Code

#### Purpose

In order to improve process control when loading or unloading the NEO Iris, many safeguards have been incorporated into the system. To indicate when it is safe to use a component, the NEO Iris uses Light Emitting Diodes (LEDs) with a consistent color code to indicate the status of specific components.

By following the consistent color code, you can continuously load and unload samples, microplates, and reagents onto the NEO Iris during operation. Continuous loading enables a high sample throughput while increasing the flexibility of the system.

Color	State	Meaning				
	Green Continuous	You can place a suitable item into this position.				
	Green Flashing	You should remove an item from this position.				
•	Orange Continuous	The system has not been requested to use this component and there are no errors.				
	Orange Flashing	A warning has been issued for this component, an error occurred or you can remove a rack.				
	Red Continuous	DO NOT place anything in this position.				
*	Red Flashing	DO NOT remove anything from this position.				

This table describes each of the LED colors within the consistent color code.

#### Color Scheme In the Plate Loading Tower

Implementation of the color scheme in the Plate Loading Tower produces the following states in the tower LEDs:



### Color Scheme In the 14-lane and 5-lane Bays

Implementation of the color scheme in the 14-lane and 5-lane bays produces the following states in the lane LEDs:



## **Barcodes**

#### Purpose

The NEO Iris uses barcodes to identify reagents, microplates, and samples in the loading bays. Using barcode technology increases the number of steps that can be automated, thus decreasing handling errors. The NEO Iris supports the use of the following barcode symbologies:

- Codabar
- Code 128
- ISBT 128 (Concatenated barcodes are not supported)
- Code 39
- Interleaved 2 of 5



**Note**: To increase data security, Immucor recommends that you use a check digit in conjunction with sample barcodes.



Limitation: Barcodes can be no longer than 18 characters in length.



**Limitation**: Barcodes must have a module size larger than 0.2mm and a bar width ratio between 2.25:1 to 3:1.

Decodability grade (grading system of A to F; A being best, F is failing) of C or better is required for consistent reading of barcodes on the instrument. This grade measures the bar width consistency throughout the barcode label. It is usually an indication of print quality of the barcode label.

The minimum number of characters in the barcode is 3.

The length of the barcode is variable but must be completely visible with a quiet zone (white space on each end of the label) of 3 mm when placed in the sample or donor rack.

The minimum height of the barcode must be 10 mm.

If barcodes have parameters outside of these specifications, barcode misreads can occur on the instrument.

Pipe characters () are not permitted as part of a sample identification.



**Limitation**: If samples have barcode identification information that is eighteen (18) characters in length and the first three (3) characters are identical to the first three (3) characters of the assay control material in assays which include plate or run controls, then the sample will be interpreted as a replicate of the control material. In this case, either the plate will fail unnecessarily (if the sample reacts differently than expected for the corresponding control), or the plate will pass but no results for that sample will be produced (if the sample reacts as expected for the corresponding control). Such a condition can also be exhibited when an assay such as crossmatch or antigen screening assay generates a set of circumstances such that the combination of donor and primary sample barcode identification information adds up to eighteen (18) characters as one of the control material barcode identifications.

#### **Reagent Barcodes**

Reagent barcodes are used to identify reagents, controls and diluents.

The following information is encoded in the reagent barcode:

- Reagent ID The type of reagent
- Lot number Batch identifier
- Expiration Date The last date that the reagent may be used
- Serial Number Unique identifier for each vial

#### **Microplate Barcodes**

Microplate barcodes are located on the short side of the plate frame. The instrument has an internal plate barcode scanner that scans the plate frame to verify the plate ID when the instrument processes the samples. For more information about the internal barcode scanner, refer to **Chapter 2 –Hardware Components**.

There are two (2) different barcode schemes used on the sides of plates, of which only one is used on a given plate frame. The barcode scheme used on a given plate frame is dependent on when a given lot of a given product is manufactured. The schemes are named *Scheme 1* and *Scheme 2*.

#### Scheme 1

The following information is sequentially encoded in the plate frame barcode for Scheme 1:

- Plate code (indicating the type of plate) linked with the lot number (batch identifier).
- Plate serial number the unique plate identifier.

Scheme 1 does not have encoded information for the plate expiration date. The plate expiration date must be entered into the software manually. Refer to **Chapter 3 – System Software Navigation** for information regarding the manual entry of plate expiration dates using the Expiry Date tab of the Plate Loading Tower dialog.

#### Scheme 2

The following information is sequentially encoded in the plate frame barcode for Scheme 2:

- Three (3) digit product identifier. For example, 008 identifies Capture-R® Select plates.
- Five (5) character expiration date (DDDYY). The DDD portion is the numeric day within the year (YY). For example, February 1<sup>st</sup> 2010 would be represented as 03210.
- Three (3) digit lot number (batch identifier).
- Five (5) digit plate serial number the unique plate identifier.

The instrument software combines the three (3) digit product identifier with the three (3) digit lot number to create the alpha-numeric lot number of the plate. The software uses a truth table to convert the product identifier to either one or two alpha character(s) that can be prefixed onto the three (3) digit lot number, the combination of which can be printed on reports. For example, SC represents Capture-R<sup>®</sup> Select plates. The printed alpha-numeric lot number of SC123 would actually be composed of 008 and 123. If a particular plate of this lot number had an expiration date of February 1<sup>st</sup> 2010 and it was the thirteenth (13<sup>th</sup>) plate manufactured, then the overall plate frame barcode would be 0080321012300013.

*Scheme 2* has encoded information for the plate expiration date. The plate expiration date does not need to be entered into the software manually. The *Expiry Date* tab of the *Plate Loading Tower* dialog will automatically populate with date information when *Scheme 2* plate frame barcodes are scanned by the instrument. Refer to **Chapter 3 – System Software Navigation** for information regarding the *Expiry Date* tab of the *Plate Loading Tower* dialog.

#### Sample Barcodes

Sample barcodes are located on sample tubes and encode the sample ID. The sample ID can be mapped to information downloaded from the host Laboratory Information System (LIS).



**Note**: For correct reading of sample barcodes in the loading bay, the sample barcodes must be positioned between 20 mm and 105 mm (0.8 and 4.1 inches) from the bottom of the tube.

## Reagents

The NEO Iris accesses all reagents, such as antisera and reagent red cells, through vials with lot numbered barcodes. The operator loads these vials into reagent racks which are slotted into the loading bay. For more information, refer to **Racks**.

The barcodes identify each vial individually so that the system can electronically monitor the fill level when the reagents are removed from the NEO Iris and then reused at a later time. When a reagent vial is empty, the software automatically moves to another bottle of the same reagent type if it is present on the instrument. When using barcoded reagents, the system automatically registers the position of the reagent so the reagent can be placed in any accessible position.



<u>Note</u>: For more information about loading reagents, refer to **Chapter 6 – Instrument Testing Operation**.

## Racks

#### Purpose

The NEO Iris uses racks for loading and unloading samples and reagents. You must place reagent vials and sample tubes in a rack and then slide the rack onto the NEO Iris. You can continuously load or unload racks from the NEO Iris during operation according to the consistent color code.

NEO Iris racks are equipped with a guide rail underneath to position them correctly in the loading bay, a handle at the back to facilitate handling, and a pin at the front that triggers a sensor to inform the system that the rack is loaded into position.

NEO Iris racks use barcodes so the system can identify them. Each rack position has its own positional barcode. On the left of the rack, next to the handle, is the overall rack barcode. The rack barcode encodes the pipetting coordinates to be used when aspirating from tubes loaded in this rack.

### Parts of the Rack

The photograph below illustrates the parts of the rack.



- A: Guide rail
- **B**: Handle
- **C**: Pin
- **D**: Positional barcodes
- E: Rack barcode

#### Types of Racks

There are two lengths of NEO Iris racks:

- Racks for the 14-lane bay
- Racks for the 5-lane bay

The following table shows every type of NEO Iris rack available in these sizes.

Bay Type	# of Positions	Rack Type	Tube or Vial Size			
14-lane bay	16	Tube: A	Sample Rack: 16–17X100 mm tube			
14-lane bay	16	Tube: B	Sample Rack: 12–13X75–100 mm tube			
14-lane bay	16	Tube: C	Sample Rack: Pediatric tube custom rack			
14-lane bay	16	Donor: D	Donor Rack: 12 mm diameter tube			
14-lane bay	12	Reagent: R	10 ml reagent vial (occupies 2 lanes)			
5-lane bay	9	Reagent: S	10 ml reagent vial (occupies 2 lanes).			
			Note: Not for use with controls.			
5-lane bay	5	Reagent: T	57 ml (43 mm diameter) vial (occupies 3 lanes)			

**Limitation**: A site visit by an Immucor representative is required to configure the Z position (downward) on the instrument for the C racks and the specific small-volume pediatric tubes in use at your site. C racks cannot be used on the instrument without this configuration. Sample probe crashes will occur without this configuration, when used in conjunction with the C racks. If differently sized small-volume pediatric tubes are subsequently used after the Z position configuration is performed using the originally designated small-volume pediatric tubes, re-configuration may be required to prevent possible sample probe crashes into the bottom of the new tubes.

## Dynamic Scheduler

The dynamic scheduler is the sequence of events that the instrument must perform to complete the requested assays. The software calculates the schedule based upon the time constraints of the individual assay steps. The system identifies or displays all resources that need to be loaded in order to complete the schedule.

The schedule is represented by the software using an axis crossing a time scale and passing through bands of color. The scheduler axis line moves in real time from left to right to demonstrate time passing as the NEO Iris operates. Each plate being run has a sequence of differently color coded bars that represent the sequence of processing steps that the plate is designed to move through for a given assay. The axis line moves through the different color bars as time progresses. The scheduler will adjust to accommodate normal processing delays.

<u> </u>	10:30 10:40 10:50	11:00 11:10	11:20 11:30	11:40 11:50	12:00 12:10	×
	R06100422					
						+
ОК						

Each band of different color bars represents one plate. The plate identification is prefixed at the left of the color band.



Each given band is spatially placed in the window to represent what the plate's starting position was in the plate loading tower. Position 1 is designated as the bottom position and position 15 is designated as the top. Therefore, for example, the plate color band for tower position 3 would be lower in the window than that for tower position 8.

Refer to **Chapter 3 – System Software Navigation** for more information regarding the schedule.



**Limitation**: Time stamps for instrument activity may not be accurate around Daylight Saving Time (DST) when a given activity spans a time period falling on both sides of the actual change of time for DST. The following recommendations are published to provide guidance on how to mitigate these time stamp inaccuracies. Allow assays to finish if they are already running during the DST change and remove the racks once processing is complete. Do not interact with the instrument (e.g. loading plates or starting assays) during the DST time change. Initialize the instrument after the DST time change is finished prior to beginning any further assays.

## Miscellaneous Information

#### **Continual Access**

The instrument is designed so that resources can be replenished and new assays started without interrupting the processing of in-progress or scheduled assays.

New samples, reagents, system liquid, and plates can be added during test processing. You can also remove completed samples, used reagent containers, plates, and liquid waste during test processing.

Instructions for replenishing samples and reagents during instrument operation are included in **Chapter 6: Instrument Testing Operation**. Instructions for replenishing system liquid during instrument operation and removing liquid waste are included in **Chapter 10: Maintaining the NEO Iris**.

#### **Expected Results**

The specific performance characteristics and the expected frequency of the possible assay results are described in the package inserts of each reagent or test well. The expected results are specific to the reagent or test wells in use. Refer to the package inserts for descriptions of expected results for reagents and test wells.

#### Biological, Electrical, Mechanical and Laser Beam Safety Precautions



**Warning**: Blood samples, liquid waste, used microplates, and consumed liquid reagent containers contain potentially biohazardous material.



**Warning**: Always wear protective gloves and clothing when handling blood samples, liquid waste, used microplates, or consumed liquid reagent containers. All blood samples, liquid waste, used microplates, and consumed liquid reagent containers must be discarded following the standard practice of the laboratory.



**Warning**: All blood products must be treated as potentially infectious. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.



**Warning**: Follow basic electrical hazard awareness to reduce the risk of injury due to prohibited electricity exposure.



**Warning**: Follow all of the necessary precautions to prevent exposure to and potential injury from instrument mechanical movement. Keep all instrument protective covers in place when operating the instrument to reduce the risk of operator injury due to instrument mechanical movement.



**Warning**: Follow all of the necessary precautions to prevent exposure to and potential injury from barcode laser scanners. Do not look directly into the laser beam of scanners or any reflections of the beam from a mirror-like surface. Exposure to the laser beam light can cause eye damage and permanent injury.
# Markings

The label is located next to the power socket on the bottom right side of the instrument. The information on the label includes, but is not limited to, the NEO Iris serial number.



# **Manufacturer Information**



Immucor Medizinische Diagnostik GmbH D-63303 Dreieich-Sprendlingen Robert-Bosch-Strasse 32 Germany

# Classification

CE

The Neo Iris Instrument is classified as "self-declared" device according to IVD Regulation (EU) 2017/746

## **Technical Support**

Please contact your local Technical Support or Immucor Technical Support International or at **+4961038056500** or via e-mail at **tech.support.int@immucor.com.** 

# **Chapter 2: Hardware Components**

### In This Chapter

This chapter describes the functions and safety aspects of the NEO Iris hardware.

For technical data about the hardware, refer to Appendix C - Hardware Technical Data.

CHAPTER 2: HARDWARE COMPONENTS	
The NEO Iris	
Cabinet	
Hood	
Loading Tower	
Transport System	
Plate Carrier	
14-lane and 5-lane Bays	
Pipetting System	
Incubator	
Washer	
Centrifuge	
Camera Reader	
Disposal Information	

# The NEO Iris

## **NEO Iris Photograph**

The photograph below shows the NEO Iris as it is viewed from the front. The instrument unit is located on top of the cabinet. The computer monitor and keyboard are attached to the right side of the cabinet.



Personal Computer (PC)

## Purpose

The NEO Iris software, installed on a stand-alone PC, controls the NEO Iris. The touch-screen PC monitor is located on the monitor table attached to the right side of the cabinet. The PC tower is located inside the cabinet.

For more information about the NEO Iris software, refer to Chapter 3 – System Software Navigation.

## How It Works

In addition to the standard PC components, the NEO Iris PC is equipped with a touch-screen monitor and a trackball (mouse).

This PC is connected to the NEO Iris, and can also be linked to an external Laboratory Information System (LIS). This enables you to download test selections (worklists) to the NEO Iris and export test results from it to the LIS.



<u>Attention</u>: Do not install additional software or add-on hardware. This will void your warranty and service contract. Addition of unapproved software may affect the performance of the NEO Iris and provide a means of introducing computer viruses.

## **Monitor Photograph**

The photograph below illustrates the touch-screen monitor and its associated parts.



## **Safety Feature**

This PC is equipped with the following safety feature:

• The computer is attached to an uninterruptible power supply (UPS) with an integral power conditioner to provide a consistent source of electricity, and to provide short periods of backup power.

# **Cabinet**

### Purpose

The cabinet is an integral part of the NEO Iris. It provides a surface on which to safely place the instrument and serves as a compartment in which to store external system components.

The section below describes the cabinet parts.

## **Cabinet Photograph**

The photograph below shows the cabinet and its parts.



- A: Common waste container
- B: System liquid containers

## **Cabinet Interior**

The following design features are incorporated into the 20 liter and 10 liter system liquid containers:

- The 20 liter system liquid container is connected to the 10 liter refill container so that when the 10 liter refill container is filled with system liquid, the liquid from that container will passively drain into the 20 liter container.
- The system is designed such that only the 10 liter refill container is manually filled (using its large cap opening), and not the 20 liter container. The 20 liter container is the direct liquid reservoir for the NEO Iris and is passively filled (indirectly) from the 10 liter refill container, using gravitational equilibration. The equilibration time for the passive draining of a full 10 liter refill container into the 20 liter container is approximately five (5) minutes.
- The 10 liter refill container can be detached from the 20 liter container (using the in-line connectors) during operation, so that it can be filled off-line, as long as there is sufficient system liquid in the 20 liter container to process the assays in progress.
- As an alternative method of filling, a commercially obtained cube of system liquid can be directly connected (using the supplied in-line connector) to the 20 liter container, instead of using the 10 liter refill container. In this instance, the commercially obtained cube of system liquid must be elevated on the drawer to at least the same level as the 10 liter refill container so that passive gravitational draining equilibration can take place. The cube must be located in the space on the drawer to the left of the common waste container and requires an extended length of tubing to connect the cube to the detachable in-line connector.

#### Cabinet





**Limitation**: Immucor requires the use of phosphate buffered (approximately 15mM) isotonic saline, pH 6.5-7.5 (PBS), on the NEO Iris system. Reactions between an antibody and its antigen may be weakened if acidic or unbuffered saline is used.

Using saline and/or deionized water in PBS preparation from sources with systems in place to control proliferation of microbes helps to reduce the chance for microbial bioburden on the system. Excessive microbial bioburden can cause degradation of system or assay performance.

## **Safety Features**

The cabinet is equipped with the following safety features:

• The 20 liter container has an integrated sensor to detect low levels of system liquid (cable circled in the photograph below).



- The base of the drawer is a sealed trough that retains any spills from the containers.
- The PC tower compartment is separated from the system liquid compartment, protecting the PC tower from splashed liquids.



<u>Warning</u>: Inadvertent operator collision with the cabinet doors or the pull out cabinet shelf can cause operator injury.

## **Cabinet Exterior**

The following design features are incorporated into the exterior of the cabinet to ensure the correct positioning of the NEO Iris.

- Two welded rings on the top of the cabinet ensure the correct positioning of the instrument and prevent it from moving during operation.
- Two rubber-tipped spacers at the back of the cabinet ensure that enough space is left for ventilation, should the cabinet be placed in front of a wall.
- Foot mountings positioned close to the wheels ensure stability of the system during operation.



<u>Attention</u>: The instrument setup process involves critical positioning and stabilization of the cabinet and centrifuge. Moving the NEO Iris could affect these aspects of instrument setup and cause centrifuge balancing errors to occur. It is recommended that the NEO Iris not be moved once installed.

# Hood

## Purpose

The hood is located on the upper front of the instrument and is used as a protective shield for the operator, to help prevent operator accidents or injuries when the instrument is performing mechanical functions. However, if required, for example in the instances of some maintenance tasks, the hinged hood can be raised to give operator access to perform necessary tasks under controlled circumstances.

## How It Works

The hood is a physical barrier that prevents operator intervention into the moving mechanical parts area of the main instrument when its components are in motion. If it is necessary to access the inside of the main instrument mechanical area, the hood can be raised up by manually pulling vertically up on the hand bar. The hood can then subsequently be lowered back into the physical barrier position by manually pulling down on the hand bar.

## **Hood Photographs**



Hood in the down protective position



Hood in the up position, for access to the internal parts of the instrument

The photograph below shows the hand bar of the hood in detail.



## Safety Features

The hood is equipped with the following safety features:

- The hood is a physical barrier that prevents operator intervention into the moving mechanical parts area of the main instrument when its components are in motion.
- The hood is linked to a sensor so that if the hood is raised under uncontrolled circumstances, such as while the mechanical parts of the main instrument are still in motion, then an audible alarm is sounded and an error message is generated which brings the instrument to an emergency stop.

# **Loading Tower**

### Purpose

You load and remove plates through the loading tower, located on the left side of the instrument and accessed from the front. You must load the plates into carrier frames first before placing them in the tower. These frames allow the transport system to move the plates and then load them into the various instrument modules required to complete the assays.

## How It Works

The loading tower is a self-contained module that enables you to load and unload up to 15 plates. The tower features a translucent door.

The plate positions are designed for easy and safe insertion of the plates. The right hand rail of each position corresponds to the guiding groove of the plate carrier, thereby prohibiting incorrect loading of the plate carriers.

Ensure also that the plates are fully pushed down into the recess of the plate carrier. When using plate strips, ensure that the strips are properly placed in their white frames and that these frames are in turn fully pushed down into the carrier frame.

A mechanical stop at the rear of each position ensures that you cannot insert the plate too far. In addition, each position has an LED indicator that shows the status of that position in accordance with the instrument color-coding system. Refer to **Consistent Color Code** in Chapter 1: Introduction to the NEO Iris for more information.

Two guidance strips attached to the inside of the door ensure correct placement of the plates in their positions when the loading tower door is being closed.



<u>Note</u>: Closing the loading tower door initiates a barcode scan of a new plate after that new plate is loaded and the LED is orange.

## Plate Loading Tower Alerts

#### Purpose

The purpose of the plate loading tower alerts is to indicate that you should not add or remove plates to or from the tower when the plate transport is accessing the tower.



<u>Warning</u>: Do not try to access the loading tower when the transport system is accessing the loading tower . You may cause a plate transport error or crash situation if you ignore the plate loading tower alerts.

#### Visual and Audible Alerts

When the plate transport is accessing the loading tower, all of the loading tower LEDs turn to continuous red.

Color	State	Meaning
Red	Continuous on all tower slot LEDs simultaneously.	DO NOT add or remove plates to or from the tower.

Additionally, during this time, if you open the loading tower front door, an audible constant tone alarm is generated that can only be silenced by the transport when it leaves the tower, or by you closing the front door before the transport leaves the tower.

It is recommended to keep the loading tower front door closed at all times, except when loading or unloading microplates.

## Loading Tower Photograph

The photograph below shows the loading tower.



## **Safety Features**

The loading tower is equipped with the following safety features:

- Integrated sensors continuously monitor the presence of plates in each position.
- LEDs indicate the status of each position at all times.

# Transport System

## Purpose

The transport system moves plates between the modules of the instrument, transporting them on frames called plate carriers. The transport system area is at the center of the instrument between the front and rear modules. In this area, the transport system can move the plates without interfering with any of the other modules. A barcode scanner is located within the transport system area, enabling automated scanning of barcodes affixed to the plates.

## How It Works

The transport system moves the plate and its plate carrier in the X, Y, and Z directions.

The transport system includes two rails that allow movement in the X-direction (left-right), one on the base of the instrument and one above the rear modules. The transport system also includes a mechanism that allows movement in the Z-direction (up and down). A Y-pusher is included in the transport system to move the plate forward and backwards, in the Y direction. Both rear and front target modules are accessible by the transport system.

During initialization of the NEO Iris, the transport system moves into home position (lower left, behind the centrifuge loading mechanism) in order to set all positioning sensors to home.

## Transport System Photograph

The photograph below illustrates the transport system and its movement.



- X: Movement
- Y: Movement
- Z: Movement

## Internal Plate Barcode Scanner

The plate barcode scanner is located behind the loading tower at the top of the instrument.



Plate barcode scanner

To scan a plate barcode, the transport system moves the plate to the barcode scanner.

## **Safety Features**

The transport system is equipped with the following safety features:

- Integrated sensors monitor all movements of the transport system and verify that the programmed steps are executed successfully. These sensors detect any obstruction in the transport system.
- The transport system only operates when the hood is closed.
- The transport system recognizes plates using a barcode system.
- The software permits manual entry of plate barcode data, with automatic logging of the operator who performs the manual entry, if plate barcodes fail to be read due to poor plate barcode position or poor printed quality of the barcode.



**Warning**: Never attempt to reach the washer area while the NEO Iris is operating. You may disrupt the instrument or injure yourself. The instrument switches off power to the motors if resistance to movement is encountered.

# Plate Carrier

## Purpose

The system transports all plates on plate carriers. This allows the instrument to process plates of different physical dimensions.

## How It Works

The plate carrier is a plastic frame containing springs and opposed mechanical stops on the inside of the frame to hold the plate firmly and correctly in place.



<u>Attention</u>: The A1 position and barcode label on the plate must face to the left (opposite the guiding groove). To avoid confusion, the upper left corner of the plate carrier is engraved with an A1. You must place the A1 position of the plate there.

The right-hand extension of the frame is comprised of a guiding groove and two holes. The guiding groove ensures correct movement and placement of the plate carrier in the modules. The two holes are engaged by pins of the transport system's Y pusher, which then pushes the plate carrier into, or pulls it out of, the different instrument modules.

## **Plate Carrier Photograph**

The photograph below shows the plate carrier and its parts.



- A: Springs
- B: Guiding groove
- C: Holes for Y pusher pins

# 14-lane and 5-lane Bays

## Purpose

The 14-lane and 5-lane bays are equipped to contain the samples and reagents required for pipetting. The bays are accessible from the front of the instrument, allowing the operator continuous access, even during instrument operation.

The 14-lane bay is in the center of the instrument and has fourteen lanes. The instrument also has a 5-lane bay, which has five lanes, on the right hand side of the instrument. The 5-lane bay is designed for use with reagents only. The 14-lane bay is designed for use with both samples and reagents.

## **How It Works**

Although both of the bays differ in size, they both work the same way. The loading bays are comprised of rack lanes, indicator LEDs, and a barcode laser scanner.

The lanes of the bays are defined by a rail that corresponds to the guide rail on the respective sample or reagent racks. Two mushroom-shaped guide buttons in front of the rail ensure that the racks are positioned correctly as they slide into the lanes. For more information about racks, refer to **Racks** in Chapter 1 – Introduction to the NEO Iris.

Each lane is equipped with an LED that indicates the lane's status in accordance with the color-coding system. For more information, refer to **Consistent Color Code** in Chapter 1 – Introduction to the NEO Iris. The indicator LED is in the first guide button of each lane.



<u>Warning</u>: Do not remove a rack when the indicator LED is flashing red, as this can damage the pipetting system and invalidates all test results on samples in the rack. Also, do not load a rack when the LED is solid red. The barcodes are not read and the reagents or samples will not be used.

The barcode laser scanner for each bay is on the right-hand side of the bay at the front and can focus on each individual lane. It reads the rack and tube barcodes as a rack slides onto the loading bay.



<u>Attention</u>: Use the correct sample and reagent racks with the appropriately sized sample tubes and reagent vials. Using an inappropriately sized tube or vial results in improper volume tracking and flagged results due to inadequate aspiration.

## 14-lane and 5-lane Bays Photograph

The photograph below shows the 14-lane and 5-lane bays.



## Safety Features

The 14-lane and 5-lane bays are equipped with the following safety features:

- Integrated sensors continuously monitor the presence of racks in the lanes and flag any rack removal or exchange.
- LEDs indicate the status of each lane at all times.
- Position barcodes on the rack prevent mix-up of sample or reagent vial barcodes and positions during rack loading.



<u>Warning</u>: Open field access is required to provide continuous access for sample/reagent loading during full system operation. Never try to access tubes or vials while their rack is still in the loading bay. Always pull their rack completely out before accessing individual positions.

Unauthorized access to the loading bay is strictly prohibited and could injure you.

# **Pipetting System**

## Purpose

The pipetting system aspirates liquids from a defined source and dispenses them in a defined destination. Both loading bays can be defined as a source location. The plate pipetting stations can be defined as both a source location and a target.

## **How It Works**

The pipetting system is comprised of the pipettors and the single probe and four probe wash stations.

## **Pipetting System Photographs**

The photographs below show the pipetting system and its parts. Under normal circumstances, when the hood is down, the pipetting system components are not visible. However, when the hood is up, the components are visible.



Four probe (left-arm)



Single probe (right-arm)

## **Pipettors**

The system has five probes mounted on two independent pipetting arms. The left arm holds four probes and is almost exclusively used to pipette samples. The right arm holds a single probe, which is primarily used for reagent pipetting. Some assays use the arms for samples and reagents.

Both of the pipettor arms are supported by a dual X-rail system. The two rails are located in the front and at the back of the main instrument body. Both pipettor arms move in the X-direction left and right across the rails.

The pipettors are fitted with steel probes. The left arm can access the entire 14-lane loading bay, but not the 5-lane loading bay. Each individual probe of the four set pipettors can aspirate up to 500  $\mu$ l. This volume can be dispensed in single or multiple deliveries. The right arm can access the entire 14-lane loading bay and the 5-lane loading bay. The single probe can aspirate up to 1000  $\mu$ l. This volume can be dispensed in single or multiple deliveries.



<u>Attention</u>: It is important to use the proper size syringe. Using the wrong size syringe can adversely affect pipetting accuracy and sensitivity. Run the Pipettor Verification Test described in **Performing the Pipettor Verification Test (PipTest)** in Chapter 10 – Maintaining the NEO Iris to confirm that the proper syringe is being used.



**Warning**: It is important to keep your hands away from the pipetting area to avoid potential injury due to the moving instrument.

All pipettors are equipped with both Liquid Level Detection (LLD) and clot detection in sample tubes, and LLD in reagent vials. The pipettors can detect volumes as low as 200  $\mu$ l in sample tubes and 1.0 ml in reagent vials. If the pipettors detect no liquid or not enough liquid, or detect a clot, the system issues a warning and, depending on assay-defined settings, waits for your correction or skips the step. The system makes an entry in the log. If low or no liquid is detected then the system issues a warning and will abort the particular plate in question and other plates scheduled to use that reagent.



<u>Limitation</u>: At least 250  $\mu$ I of packed red blood cells need to be present in a sample tube to ensure that the probe picks up red blood cells and not plasma (only for those assays that require red blood cells).

At least 500 µl of plasma or serum needs to be present in a sample tube to ensure that the probe picks up plasma or serum.



<u>Attention</u>: If the system detects a clot, it is important to check that the probe is not blocked, as this inhibits further pipetting and leads to operational failures. For more information, refer to **Chapter 6 – Instrument Testing Operation**.

During instrument initialization, both the pipettor arms and the pipettor pumps undergo a full motion check, during which they move into their home positions to set all positioning sensors to zero.

## **De-bubbler Module**

The de-bubbler module is an integral part of the pipetting system.



De-bubbler module

The following design features are incorporated into the de-bubbler module.

- The de-bubbler module is located inside the NEO Iris, on the rear right-hand back corner, to the right of the dilutor syringes, as the instrument is viewed from the front.
- The de-bubbler module is a vertical clear plastic cylinder with tubing entering and exiting the cylinder from above and below.
- The de-bubbler is designed such that system liquid is pulled into the chamber by a pump that draws from the tubing at the top. Any bubbles that enter the chamber rise and are pulled out via the top tubing and are returned to the 20L system liquid container. The liquid at the bottom of the chamber remains free of bubbles. It is this liquid that is drawn from the bottom of the chamber that is used in the pipetting system.
- No PBS is wasted by the de-bubbler module during assay processing.

## **Probe Wash Stations**

To prevent contamination due to carry-over between samples, the steel probes of the sample pipettors are washed in the wash stations after every dispense step.

You can also use the wash station to prime the liquid system prior to its first use.





Four Probe Wash Station

Single Probe Wash Station

The wash stations consist of an elevated plastic cup containing one or four rinse wells. Washing is a two-step process, as follows:

- The pipettors move to the waste position over the wash station cup where the remaining contents of the probes are emptied by flushing system liquid through and out of the probe.
- The pipettors move over and into the rinse wells and system liquid is pumped through the pipettor tubing. This action washes the inside and the lower outside end of the steel probes simultaneously.



<u>Attention</u>: The instrument cannot differentiate between water and system fluid. If deionized water is mistakenly used as the system fluid, the test results are invalid.

<u>Attention</u>: For the probe wash to work properly, the settings of the instrument, with regards to proper positioning of the probes inside the rinse wells, are crucial. For details, refer to **Checking Pipettor Reference** in Chapter 10 – Maintaining the NEO Iris.

After sampling red blood cells from the bottom of sample tubes, the system performs a special external rinse in the external rinse positions. The probes move to the bottom of the external rinse wells and the active wash pump rinses the outside of the probes by pumping system liquid. The probes move slowly out of the external rinse wells to achieve complete rinsing of the outside of the probes.

The system actively pumps waste liquid from the wash stations to the common waste container in the cabinet.

## Safety Features

The pipetting system is equipped with the following safety features:

- Integrated sensors monitor the position of the pipettors as well as pipettor movement. If something hinders pipettor movement, the system displays an error message and movement stops.
- The pipettor arms are supported by a dual X-rail system that provides stability for the moving components.
- The probes are spring loaded to absorb minor bumps so that damage is minimized, should the probes be driven with force into an unexpected obstacle.
- Liquid level detection, clot detection, and aspiration tracking ensure accurate aspiration of liquids or adequate flagging of results.
- Syringes are visible to allow for inspection.
- Software controls limit collision of both pipetting arms.
- If the hood is opened while the instrument is processing, any pipettor arm movement stops and the instrument needs to be re-initialized.
- You can verify proper pipetting by observing correct reciprocal forward and reverse ABO hemagglutination reactions and also detecting red blood cells in test wells where they are not supposed to be delivered. Such inappropriate delivery of red blood cells could be due to, for example, when the plasma volume is too low and the probe picks up red blood cells instead of the intended plasma.
- Software controls monitor the addition of red blood cells, plasma (or serum), and reagents to test wells.
- After sampling red blood cells from the bottom of sample tubes, the system performs a special external rinse in the external rinse positions. The system also rinses the inside of the probe to remove plasma, red blood cells, and reagents in order to prevent contamination.
- The de-bubbler module removes any bubbles from the system liquid used in the pipetting system.
- Software monitors processing and post-processing allowable times for time sensitive assay steps and will invalidate results if exceeded.
- A filter in the liquid supply prevents the system from pipetting particulate material. Clogged probes result in high-pressure errors.

# Incubator

## Purpose

The incubator provides the appropriate environment for the incubation steps of the assays performed on the instrument. The incubator is located to the right of the reader at the back of the instrument.

## How It Works

The incubator comprises 15 individual positions. These are divided into two different temperature zones—a room temperature zone (6 positions) and a 37°C zone (8 positions)—and a neutral insulator zone that separates the two different temperature zones.

The thermodynamics of the 37°C incubator bays require an actual temperature of 39°C to yield a test well temperature of 37°C. Each heating plate incorporates its own software controlled temperature sensor. The position between the temperature zones is intended as an insulator and cannot be used for incubation. This insulator is an active peltier element that inhibits heat transfer from the warmer to the cooler zone.

The incubation positions are heated from above and below by resistance foils integrated into the aluminum plates that separate the positions.

A flap, to retain heat, covers the opening of each position. Each flap opens individually to allow the transport system to slide in or remove a plate. Plate carriers are held securely inside the incubator by magnets at the back of each position.



**Limitation**: Laboratory ambient temperature and humidity affects the room temperature incubation bays, and an elevated ambient room temperature may disrupt assays that need to be incubated at specified temperature ranges, as published in the relevant package insert.

## **Incubator Photograph**

The photograph below shows the incubator.



## **Safety Features**

The incubator is equipped with the following safety features:

- Integrated sensors monitor the temperature and the presence of plates in each position.
- Sensors monitor the open/closed status of the front flaps.
- The heating power of the foils is limited so that even after a severe error leading to complete loss of temperature control, the system can only reach a maximum temperature of 70°C.
- Instrument software monitors and controls incubation time.
- The temperature distribution over a plate during incubation is less than 1°C from average. Results are invalidated if temperatures exceed the desired range.
- The instrument software does not permit the start of an assay that requires ambient temperature incubation if the incubator is too hot or too cold.

# Washer

## Purpose

The washer performs all plate-washing cycles required for an assay protocol.

## Washer Module

The washer, located in the back of the instrument, comprises the wash head assembly, the pumps, the plate loading frame, and the overflow and prime troughs.



The transport system slides plates in plate carriers into the washer. Two magnets at the back of the washer ensure that the plates are correctly and securely positioned.

The wash head assembly, which moves left to right above the plate, aspirates liquid from the plate wells, and dispenses clean system liquid into them as required by the assay protocol. The wash head assembly comprises eight pairs of aspirating/dispensing needles. The aspirating needles are longer than the dispensing needles. The wash head assembly is lowered into the plate wells for each aspiration or dispensing step.

The priming trough is located on the left-hand side. The wash head assembly dispenses into this trough during washer priming. Aspiration of the priming liquid runs continuously through the aspiration function of the wash head assembly is not in use, it is positioned in the priming trough. The priming trough includes a Liquid Level Sensor (LLS) board that verifies liquid is at the expected level during aspiration and dispense checks. An alarm will be generated if the LLS detects an abnormality.

If a spill occurs because of an error, the overflow trough catches the spilled liquid. The overflow trough is equipped with a sensor mat that issues an alarm and stops the washer should the liquid fall into the trough during processing.

## Safety Features

The washer is equipped with the following safety features:

- An integrated sensor monitors the pump for correct operation.
- A priming trough that is fitted with a Liquid Level Sensor (LLS) to ensure that spilled liquids are detected. The LLS is involved in wash verification and the detection of manifold failures.



- An overflow trough that catches spilled liquid.
- A liquid overflow detection mat is located below the plate washing area so that manifold fluid leaks over the washing area can be detected. A software error is generated if fluid is detected.
- Closed containers house potentially bio-hazardous waste materials. All connectors that could come in contact with bio-hazardous waste are fitted with check valves to prevent liquids from dripping upon disconnection.
- A priming step to flush bubbles from the lines prior to each plate wash.
- A manifold aspiration and dispense verification check is performed before and after every plate wash to make sure each well is washed and aspirated.
- Monitoring of system liquid levels that alerts you when the volume is low.
- The system monitors waste liquid level and alerts you when the volume is high and the container needs to be emptied.
- The instrument will not use the system liquid if the level sensing device is not connected.
- Monitoring of plate positions into and out of the washer. You are alerted of any jams.

# Centrifuge

### Purpose

The centrifuge module is located on the left-hand side of the NEO Iris as a detachable unit and performs two functions:

- The centrifuge spins the plates with g-forces up to 1200 xg.
- The centrifuge shakes the contents of the plate wells to resuspend sedimented red blood cells.

The centrifuge-loading unit to the left of the loading tower is a functional part of the centrifuge that transfers the plate from the main instrument into the centrifuge.

## How It Works

**Centrifugation** - Plates can have different weights depending on the assay performed and whether the plate is fully or partially used. The centrifuge must be balanced at all times. The system achieves the balance by constantly monitoring the vibrations caused by the centrifugation. If the vibrations exceed a certain threshold, the system performs a balancing step during which it adjusts a counterweight to minimize vibration level while the centrifuge is spinning at constant speed. The system then further accelerates to the assay-defined g-force. The rotor is the arm in the centrifuge that spins the plate.

Shaking – The centrifuge module can also shake plates in order to mix the well contents.

## Centrifuge Photograph

The photograph below shows the centrifuge.



# **Decoupling of the Centrifuge**

Centrifuge decoupling describes the physical separation of the centrifuge module from the main instrument. The decoupling of the centrifuge from the main instrument is done to alleviate centrifuge balancing problems that subsequently generate periodic balancing timeout errors. These timeout errors stop processing.

Decoupling is visible to you by the presence of

• a metal base plate beneath the centrifuge (circled in yellow in the photograph below), and



• black rubber feet in all four (4) bottom corners of the centrifuge.





**Note**: The rubber feet provide a flexible mounting for the centrifuge.

The use of decoupling helps to provide

- better vibration adjustment by the NEO, and
- reduction of vibration transfer to the main instrument, thereby reducing such vibrations from being transferred to the pipetting system and other modules of the main instrument.



**Note**: If you encounter balancing timeout errors, you should call Technical Support to report them regardless of the fact that your centrifuge has been decoupled.

The following warning applies to your interactions with the centrifuge that has been decoupled from the main instrument.



**Warning**: You must not bump, knock, rest up against or otherwise come into physical contact with the centrifuge module because you may cause balancing errors or centrifuge loading and unloading issues. By causing these errors or issues, you could therefore interfere with assay processing. No items should be stored on top of the centrifuge cover including, but not limited to, sample or reagent racks, books, documents or other laboratory consumables and supplies.

## Safety Features

The centrifuge is equipped with the following safety features:

- The centrifuge is entirely enclosed within a containment barrier during operation.
- The centrifuge rotor housing withstands complete rotor disintegration at full speed with only minor dents. No parts will leave the housing.
- Integrated vibration sensors continuously monitor the centrifuge. Excessive vibrations activate the emergency brake of the centrifuge.
- Sensors monitor the centrifuge loading door and service hatch. The centrifuge cannot spin unless both doors are closed.
- The software monitors centrifuge speed. It flags too high and too low speed conditions.
- The centrifuge loading unit (CLU) has a plexiglass window in the top of the cover so that any CLU problems can be viewed by the operator without the need to remove the CLU cover.



• The CLU cover is fixed by a single rounded Allen head screw and an associated washer, in addition to VELCRO<sup>®</sup> fasteners, so that it is only removable by using a tool (2.5 mm Allen key). This is to prevent operator access to the inside of the centrifuge in the event of complete loss of software control leading to the centrifuge spinning and the loading door falling open. The screw (with the washer) is located on the lower middle of the front of the cover.



- The centrifuge is located on top of the cabinet with small foot risers to incorporate additional stability.
- The software will alert the user if there is an error loading or unloading a plate into the centrifuge.

# **Camera Reader**

## Purpose

The reader module is located in the rear of the instrument and takes an image of the plate.

## How It Works

The reader uses Complementary Metal-Oxide Semiconductor (CMOS) cameras to take an image of the plate from underneath. The system calculates a value for each well based on a multi feature image analysis.



**Limitation**: The NEO Iris must be switched on at least 30 minutes prior to the first plate read to allow the reader lamp to warm up. Reading of plates prior to completion of this warming period can cause incorrect negative reading of weakly positive reactions.

## Camera Reader Photograph

The photograph below shows the camera reader position with the door closed.



## **Safety Features**

The camera reader is equipped with the following safety features:

- Potential air bubbles and foaming in the wells (where the camera takes an image) are minimized by the pipetting system, by use of the de-bubbler module and optimized aspiration/dispense procedures. Air bubbles in the test wells are located on the perimeter and do not interfere with reading of results.
- If the reader is unable to take an image the plate will automatically be aborted.
- The software recognizes incorrect plate placement and will not take an image if it is in the incorrect orientation.

# **Disposal Information**

## Liquid Waste

All liquid waste is collected in the common waste container located in the cabinet.



<u>Warning</u>: The liquid waste is potentially biohazardous material. Always wear protective gloves and clothing when handling the liquid waste. If any liquid waste is spilled, clean it up immediately following the standard practice of the laboratory.

Empty the common waste container using the procedure described in **Chapter 10 – Maintaining the NEO Iris**. Dispose of the contents according to the laboratory guidelines.

## Solid Waste

The only solid waste produced by the NEO Iris is used plates and consumed liquid reagent vials.



<u>Warning</u>: Used plates and consumed liquid reagent vials contain potentially biohazardous material. Always wear protective gloves and clothing when handling used plates. If any liquid from a plate is spilled, clean it up immediately following the standard practice of the laboratory.

Plates are returned to the loading tower after they have been processed. Remove the plates from the loading tower and dispose of them according to the laboratory guidelines.



<u>Note</u>: Only remove plates from the loading tower when the indicator LED for that position is orange or flashing green.

Empty reagent vials can be disposed of after removing the reagent rack from the instrument.



<u>Note</u>: Only remove racks from the loading bay when the indicator LED for that position is orange or flashing green.

# Chapter 3: System Software Navigation

## In This Chapter

This chapter describes the NEO Iris software.

CHAPTER 3:	SYSTEM SOFTWARE NAVIGATION	3-1
Navigation	۱	3-2
Main Men	u Bar	3-8
Machine M	1onitor	3-59
Status Bar		

# Navigation

## Screen Layout

### Main Screen Components

The instrument's graphical user interface, the Main Screen, reflects both the structure of the instrument and the pattern of the workflow. This makes it easier to navigate through the software as well as providing important information when relevant.



The screen is divided into three areas:

Main Menu Bar

Machine Monitor

Status Bar

#### Main Menu Bar



The Main Menu Bar at the top of the screen provides access to system dialogs. You use the system dialogs to program system settings, perform maintenance, start an assay run, log in to and out of the system, and exit the system.

#### **Machine Monitor**



The Machine Monitor, below the Main Menu Bar, is the primary display area.

The Machine Monitor is a graphical top view of the instrument that provides an overview of the status of the modules and access to the module dialogs. You use the module dialogs to view and edit information about the module. For complete information about the Machine Monitor and the module dialogs, refer to Machine Monitor later in this chapter.

#### **Status Bar**



The Status Bar at the bottom of the screen is always visible, and provides an overview of the system status. The Status Bar includes, from left to right, a Log List button, the NEO Iris Emergency Stop button, some wording that indicates the overall status of the NEO Iris and three system status indicators.

The Log List button on the left side of the status bar displays a list of the last actions performed by the system. Press the Log List button to display the Log List dialog, which displays a full log list.

For more information about the Log List, the Emergency Stop button, the wording that indicates the overall status of the NEO Iris and the three system status indicators, refer to **Using the Status Bar** later in this chapter.

#### **Confirmation Dialogs**

Confirmation dialogs are safety features in the operating software. The system displays these dialogs after you request certain actions.

2	Do you really want to initialize?	
~ L		
Γ		
	Yes	No

The software prompts you to confirm that you wish to continue with the action. You must either click **Yes** if you want to continue or click **No** if you do not want to continue.

### **Common Buttons**

The following buttons are standard command buttons used throughout the software. They always have the same function.

Button	Description	
ОК	Close the dialog and save any changes made.	
Cancel	Close the dialog without saving any changes.	
Done	Close the dialog.	
Print	Print information from the dialog.	

## **Navigation Options**

#### Keyboard and Trackball

The PC uses a standard keyboard and trackball mouse. You can use the keyboard to type alphanumeric characters into the field and to navigate through the software.

To navigate through the software using the keyboard, press the **Tab** key on the keyboard to move from one element in the dialog to the next. Press **Enter** to execute the action of the element currently selected.

Selected?	Button	Text Box
Not selected	ОК	
Selected	ΟΚ	

You can press the **Print Screen** key on the keyboard to print a paper copy of whatever area of the instrument software you are currently accessing on the computer monitor.

#### Touch Screen

The touch screen enables you to interact directly with the software. Most elements that you can access using the keyboard and mouse can also be accessed using the touch screen.

To use the touch screen, simply touch the element or button on the screen that you want to use, such as a menu bar button, list box, text box, or action button.



**Note**: When touching the screen, be sure to tap or press gently, and do not press excessively. Excessive or prolonged cumulative screen pressure can misalign the screen's responsiveness to touch or cause other issues such as software freezes.

#### When to Use

Throughout this manual, when you are instructed to click or press an element or button, you can do one of the following:

Touch the element or button directly on the screen.

Use the trackball to place the pointer over the item and click the left trackball button.

Use the **Tab** key to move through the dialog box to the element needed and press **Enter**.

For example, to display the *Login* dialog, click (or touch) the **Login** button on the Main Menu Bar.
# Main Menu Bar

## About the Main Menu Bar

## Main Menu Bar icons

The table below describes the icons in the Main Menu Bar.

lcon	Name	Description
5	Login	Allows you to log in to the system.
	Initialize Instrument	Allows you to initialize the instrument.
\$	Start Run Assistant	Provides an intuitive guide throughout the procedures necessary to start an assay run.
	Worklist Editor	Allows you to view, edit and request tests for samples in the software.
	Plate Status List and Graphical Schedule	Displays data about all plates that have been processed. You can also search for, view and delete plate data.
Results	Result Viewer	Allows you to view test results.
	Maintenance	Allows you to view and activate automated maintenance checks.
X	Utilities	Allows you to view event logs and statistics, as well as archive and print reports.
*	Instrument Settings	Allows you to program the system settings in order to configure the instrument to the individual needs of the laboratory.
<b></b>	Version Information/Operator Manual/Help	Displays the electronic operator manual on the instrument computer. Also displays a list of the modules and software used on the instrument.
<b>-</b> 3]	Shutdown	Guides you through the procedure that you must perform before the system is shut down.

# Logging In

The Login dialog allows you to log in to the system. Click the **Login** button on the Main Menu Bar to log in.



The system displays the Login dialog.

NE	) Iris - Login
Γ	Please enter your name and password
	Name :
	Password :
	Start Action Restore Databases from Backup
	Use DMS Archive
	0K Shutdown

Once you are logged in, you can click the Login button again to log out. The system also displays the Login dialog during system startup.

The **Login** dialog contains the Name and Password fields. As soon as you type at least one character in both fields, the dialog activates the OK button.

The **Shutdown** button is only accessible at first login. This allows you to shut down the computer if start up was accidental.



**Note**: The system can be configured to perform an automatic Logoff (Option available in Instrument Control Software version 1.8 and higher). The software will activate an **AutoLogoff** locking screen after a certain period of inactivity (in minutes). The duration of the inactivity period prior to Auto Logoff is configured by Immucor personnel.

A	utoLogoff	
	Please enter you	r name and password
	Name :	
	Password :	
		OV

The **AutoLogoff** dialog contains the Name and Password fields. As soon as you type at least one character in both fields, the dialog activates the OK button.

Click the **OK** button to unlock the system.

For information about how to log in to the instrument, refer to **Logging In, Chapter 5 – Instrument Start-Up**.

For information about setting user names, passwords, and access rights, refer to **Assigning Passwords** and User Access Rights in Chapter 4 - Security.

## Initializing the Instrument

The **Initialize Instrument** button allows you to initialize the NEO Iris. Initialization is used to perform critical equipment verifications and resets prior to allowing the NEO Iris to be used for the necessary maintenance tasks and subsequent assay processing. Initialization resets the whole system, primes all tubing, returns all modules to their home positions, and moves any plates on the NEO Iris back to the plate loading tower.

Press the Initialize Instrument button to initialize the NEO Iris.



The system displays the *Confirm Initialization* dialog for you to confirm initialization.

Confirm 1	Initialization	
2	Do you rea	ally want to initialize?
[	Yes	No

The *Confirm Initialization* dialog is closed by pressing the **No** button. The *Instrument Initialization* dialog is displayed by pressing the **Yes** button of the *Confirm Initialization* dialog. The *Instrument Initialization* dialog sequentially lists all individual initialization activities in the *Action* list as they are performed by the NEO Iris.

Action		
initializing COP initializing devices Pipettor ok Transport Incubator1-{003F} - ok Washer1 Centrifuge BarcodeReader - ok Incubator1-{7F80} - ok Tower - ok Sample Loading Bay IOModuleAdapter - ok DataReduction - ok		4
	Cancel	

Once initialization is complete, the *Instrument Initialization* dialog is no longer displayed. Initialization can be stopped prior to its completion by pressing the **Cancel** button of the *Instrument Initialization* dialog.

The system automatically initializes each time you turn on the NEO Iris. However, it may be necessary to request initialization manually, for example, after leaving the instrument in standby for a long period, or as a last method to recover from a fatal error situation that you cannot otherwise correct.



**Note**: During initialization, the system deletes all reagent and sample rack information. You must reload the racks before starting the next run.

## **Start Run Assistant**

#### About the Start Run Assistant

The *Start Run Assistant* dialog provides an intuitive guide through the procedures necessary to start an assay run. Click the **Start Run Assistant** button.



The system displays the Start Run Assistant dialog.

	Start Run Assistant	
. If specimen tubes have not yet been loaded please g	o to "Load Samples" .	Load Samples
2. Query the host system for test orders.		Download Requests
<ol> <li>Complete loading (micro plates, reagents, controls, dis</li> </ol>	posable tips, and wash buffers).	Load Resources
Cancel		

The Start Run Assistant dialog presents you with three preparation stages:

Load samples.

Download requests from the host Laboratory Information System (LIS).

Load Resources.

Refer to Chapter 6 - Instrument Testing Operation for more information about the Start Run Assistant.

#### **Resource Overview Window**

The *Resource Overview* window provides an overview of the consumables required to run all the assays requested for the samples loaded on the instrument. Click the **Load Resources...** button on the *Start Run Assistant* dialog to display the *Resource Overview* window.

Resources	Assay Name	Samples/Strips	Plates	Reagents	Controls	Donors	Washbuffer	Pipettor	Incubator
	DAT	1/1							
	Pool_Cell	1/1							
0	ReflexABO	1/1	<b></b>	0	÷	-	÷		<b></b>
0	ReflexFWD	1/1	0	0	÷	÷	÷	<b>V</b>	÷
Star	i:	Cancel							

#### How It Works

The system groups all test orders for the loaded samples onto plates. Each line in the dialog corresponds to one plate. Based on the Samples/Strips column, you can decide to start the plate or not.

To select a plate to be run, click in the corresponding line. The line highlights blue, and the system performs a resource check for all consumables needed to run that plate. If the system displays a red exclamation mark after selecting an assay, you can check the corresponding resource by clicking the button above each column. If you can use the resource loaded on the system, the system displays a green checkmark  $\blacksquare$ .

During scheduling, or if the resource is not used for the requested assay, the system displays a gray division sign in the column below the named resource.



#### **Resource Overview Window Columns**

The *Resource Overview* window displays a table with ten columns. When all Resources fields are either blank or have a green check mark displayed, you can click the Start button to start the run, otherwise the **Start** button is unavailable. The system generates a processing schedule (after start is selected) and a shopping list (if an item is deficient). Both are available to you.

The table below describes each of the columns in the *Resource Overview* window:

Column	Description
Resources	A summary of the resource check. If all required resources are present, the system displays a green checkmark 🗹.
Resources	If resources are missing or user maintenance is due, the system displays a red exclamation mark.
Assay Name	The name of the assay to be run.
Samples/ Strips	The number of samples and strips assigned to that plate.
	If the assay is sample driven, the system displays each sample request in this list. If there is no plate present, the system displays a red exclamation mark.
Plates	If a plate is present, the system displays a green checkmark 🗹
	If the assay is plate driven, every plate loaded in the loading tower generates a request for that plate. Clicking the <b>Plates</b> button at the top opens the <b>Loading Plate</b> dialog.
Reagents	If all reagents necessary for the assay are loaded, the system displays a green checkmark $\square$ .
	Clicking the button at the top accesses the 5-lane reagent loading bay.
	If all controls are loaded, the system displays a green check mark.
Controls	
	Clicking the <b>Controls</b> button accesses the 14-lane loading bay.

Column	Description
Donors	In all assays except the crossmatch and antigen screening assays, the system displays a gray division sign. If all requested donors are present for an assay, the system displays a green check mark.
	Clicking the <b>Donors</b> button accesses the 14-lane loading bay.
Machbuffor	If the system displays a red exclamation mark 🔟, you must check for the presence or absence of system liquid.
Washbuller	You must also click the <b>Initial Prime</b> button on the <b>Wash Buffers</b> dialog for the specified buffer container.
Pipettor	If the system displays a red exclamation mark, the system liquid is empty or the waste is full.
	In this case, refill the system liquid or empty the waste. Clicking the <b>Pipettor</b> button displays a more detailed description of what actually caused the error.
Incubator	If the system displays a red exclamation mark, the temperature in the incubator slots is out of range for the temperature zone required in this assay.
	Clicking the Incubator button accesses the Incubator dialog.

## Work List Editor

#### About the Work List Editor

The *Work List* dialog allows you to manually view, edit, and request tests for samples in the software. Click the **Work List** button on the Main Menu Bar.



The system displays the Work List dialog.

Reference	Sample ID	Assay	Donation ID	Units Requested	Requested by	Requested on	Transferred	Priority	First name	Last
0001	243168070301090672	Ab_ID			yolanda	2009/08/07	1			
0002	243168070301090672	Ab_ID			yolanda	2009/08/07	1			
0003	243168070301090672	Ab_ID			yolanda	2009/08/07	1			
0004	R10623	ReflexABO			yolanda	2009/08/07	$\checkmark$			
0005	L.	Weak_D_F			yolanda	2009/08/07	1			
0006	Eisenschenk	Weak_D_F			yolanda	2009/08/07	1			
0007	R10623	ReflexFWD			yolanda	2009/08/07	1			
•										•
			Add	X Delete	Refresh	Donation	ns <u>E</u>	<u>Print</u>	ок	

Column headers are displayed that describe various facets of the worklist entries. The column headers are shown below.

Reference	Sample ID	Assay	Donation ID	Units Requested	Requested by	Requested on
Transferre	d Priority	First name	Last name	Gender Date o	f Birth	

Patient demographic details will only be displayed under the relevant columns for a given sample if the patient already has one or more sample results stored on the computer hard drive and those demographic details have been manually entered on the *Test Overview* tab display of the sample result *Test Details* dialog for those previous sample(s). The **Gender** drop-down list allows for selection of *Male* or *Female*. The **Date of Birth** drop-down list allows for the selection of a date from a navigable calendar. The *Last name* of the patient name should be entered into the first field of **Patient Name** and the *First name* should be entered into the second field of **Patient Name**. Additional fields of **Patient Name** allow for the entry of middle names or other naming conventions.

Test Querview	Descents   Euro	t and Distance		
results	Reagents   Ever	ic Log   Place vie	ews	
Sample Details				
Sample ID : LQ13608	3		Gender :	
Patient Name :				
Patient ID :			Date of Birth : 2010/0	3/01

The window of the *Work List* dialog displays the work list entries. The order of these entries can be changed if one of the header descriptions is preferred as a sorting mechanism. For example, by pressing the **Sample ID** header, the order of all of the entries in the window will be rearranged according to alpha-numeric sample order.

The activity buttons are located at the bottom of the *Work List* dialog.



**Note**: The *OK* button is inactive when the *Refresh* button is pressed, and it remains inactive until the *Refresh* button task is completed.

#### **Work List Editor Buttons**

The table below describes the buttons in the Work List Editor.

Button	Option	Description
STAT	STAT	Facilitates the ability to assign status of STAT to work list entries.A small version of the icon will be displayed under the <i>Priority</i> column if this status is applied to a specific sample or samples.Image: Mote: STAT priority cannot be removed using this button.
Add	Add	Display the Worklist–Add Items dialog.
<b>X</b> Delete	Delete	Delete a selected item. The system prompts you for confirmation. The default button in this dialog is the <b>No</b> button. Pressing the <b>Enter</b> or <b>Esc</b> key while in the dialog results in nothing being deleted. When no items are selected, the <b>Delete</b> button is unavailable. Pressing the <b>Delete</b> key on the keyboard has the same effect as pressing the <b>Delete</b> button.

Refresh	Refresh	Refresh the worklist. This has no effect on the automatic, timer operated refresh operations or their frequency. The <b>Refresh</b> button is always available.
Donations	Donations	Add or delete donor IDs for crossmatch tests. This button is only available when you have selected a crossmatch test master entry, which specifies the sample ID and the number of donations. When you click the <b>Donations</b> button, the system displays the <i>Worklist–Add Items</i> dialog.
Print	Print	Print the contents of the worklist. The system displays the <i>Print</i> dialog. If you have selected lines, the system prints only those lines. If you have not selected any lines, the system prints the entire worklist. The <b>Print</b> button is always enabled.
ОК	ОК	Method to exit dialog back to the Machine Monitor.

## Worklist-Add Items

The system displays the *Worklist–Add Items* dialog when you click the **Add** button from the worklist. Use this dialog to specify new worklist entries.

₩orklist - Add It	ems					×
Sample ID :	LS061502	Dona	ition Count 🛛 🛨	Donation ID	:	
Tests Requested				-3). [		
FWDABO_AB	IgG_XM	mAgScr37	mAgScrRT			
pAgScrAHG	Pool_Cell	ReflexABO	ReflexFWD			
Rev_ABO	RfxABO_AB	RfxFWD_AB	Weak_D			
Weak_D_F	Wk_D_AB	Wk_D_F_AB			7 -24-	
				Prev Next	Remove Donation	STAT
Custom Profiles -	2 3	4	5	Line Canala	En la Fuit	×

In the case of crossmatch and other similar assays, the system disables the Sample ID, test selection buttons, all **Custom Profiles** buttons, and the **Next Sample** button in the *Worklist–Add Items* dialog. Only Donation Count, Donation ID, donation list box, **Remove Donation**, **Save and Exit**, and **Cancel** are available. The title of the dialog becomes *Worklist–Add Donations* instead of the default *Worklist–Add Items*. The Donation ID is the first field with focus.

### Worklist-Add Items fields and buttons

The table below describes the fields and buttons in the *Worklist – Add Items* screen.

Field	Description	
Sample ID	Use the Sample ID field to identify the sample. The system automatically displays the cursor in this field. It contains the sample ID of the selected line from the <i>Worklist</i> dialog if any lines were selected when you clicked the <b>Add</b> button from the worklist. If no lines were selected, the input line for Sample ID is blank.	
Donation Count	Use the Donation Count field to specify how many donors can be crossmatched. The default for this box is blank. The field has an up-down control with a range of 0 to 99. The control also works with the arrow cursor keys. This field is disabled if no crossmatch test is selected.	

Field	Description	
Donation ID	Use the Donation ID field for crossmatch testing by entering the donor unit identification numbers. Beneath the entry is a list box that contains the donation list. The list box is a multiple-selection type list box. This field is disabled if no crossmatch test is selected.	
	Press the Enter key from the Donation ID text box to perform the following sequence:	
	1. The system adds the item in the Donation ID field to the multiple-selection list box below it. If the entry already exists in the list box, the system does not add it again, but highlights the existing entry.	
Enter	2. The system highlights the new item in the list box.	
	3. The system clears the Donation ID field.	
	4. The system places the text cursor in the Donation ID field.	
	5. If adding a donor causes the number of donors to exceed the donor count, the	
	system displays a warning dialog box.	
Remove Donation	Click the <b>Remove Donation</b> button to delete the selected entries. If the list box is empty, or no list item is selected, the <b>Remove Donation</b> button is unavailable.	
STAT	Use this button to assign the status of STAT to the added entry.	
Test Buttons	The <i>Worklist – Add Items</i> dialog contains a special button for each test. Each button contains the short name for the test and an indicator light. The light is red for unselected tests and green for selected tests. Additionally, the button is in a pressed state when selected and a raised state when not selected. All test buttons are initially red when the system opens the dialog. Pressing a button toggles it between selected and not selected.	
Next/Prev	If there are more tests than can fit in the dialog, the <b>Next</b> and <b>Prev</b> buttons enable you to scroll through the complete set of available tests. The <b>Next</b> button is disabled when there are no more tests to display. The <b>Prev</b> button is disabled when there are no preceding tests to display.	
Custom Profiles	The dialog contains <b>Custom Profiles</b> shortcut buttons that control which assays are selected. Click one of these buttons to select a particular profile of assays. When yo click a custom profile button, the system selects every assay associated with that custom profile and excludes all other assays.	

Field	Description	
	Click the <b>Next Sample</b> button in the Worklist–Add Item dialog to perform the following actions:	
	1. The system saves the current information to the worklist.	
	2. The system clears the Sample ID.	
	3. The system places the text cursor in the Sample ID field.	
Next Sample	4. The system clears the Donor Count.	
	5. The system clears the Donor ID.	
	6. The system clears the Donations list box.	
	7. The system deselects all tests (red).	
	The Next Sample button is unavailable whenever the Sample ID is blank or no tests	
	are selected.	
Save & Exit	Click the <b>Save &amp; Exit</b> button to close the dialog, save the current changes, and add the request to the worklist. The <b>Save &amp; Exit</b> button is unavailable if the Sample ID is blank or if no check boxes are selected.	
Cancel	Click the <b>Cancel</b> button to close the dialog without saving the current changes. Pressing the <b>Esc</b> key has the same effect as pressing the <b>Cancel</b> button.	

## Plate List

## About the Plate List

The *Plate List* dialog displays data about all plates that have been processed. You can also use this dialog to search for, view, and delete plate data. Click the **Plate List** button.

The system displays the *Plate List* dialog.

Plate ID	Assay Code	Status	finished at 💌	last Event
FC350	Flow350	finished	2009-08-11, 15:00:54	plate moved to PlateTower (2009-08-11, 15:00:53)
PC0000009	PCheck	finished	2009-08-11, 14:41:52	Flag Evaluation: Passed (2009-08-11, 14:41:52)
PC0000009	PCheck	finished	2009-08-11, 14:35:05	Flag Evaluation: Failed (2009-08-11, 14:35:05)
gfcghcb	PrimeAir	finished	2009-08-11, 13:37:51	plate moved to PlateTower (2009-08-11, 13:37:49)
R04802092	3_Cell	finished	2009-08-11, 12:10:00	plate moved to PlateTower (2009-08-11, 12:09:59)
R04802092	3_Cell	finished	2009-08-11, 11:34:23	plate moved to PlateTower (2009-08-11, 11:34:22)
R04802092	QC3_Cell	finished	2009-08-10, 15:56:31	Evaluation status set to ok. (2009-08-10, 15:56:37)
R04802092	QC3_Cell	aborted	2009-08-10, 15:21:58	plate moved to PlateTower (2009-08-10, 15:21:57)
RESV	ResVol	finished	2009-08-10, 14:58:23	plate moved to PlateTower (2009-08-10, 14:58:22)
RESV	ResVol	finished	2009-08-10, 14:53:43	plate moved to PlateTower (2009-08-10, 14:53:41)

The *Plate List* dialog is divided into three areas:

Information area at the top

Plate status list in the middle

The buttons at the bottom

### **Information Area**

The information area at the top of the dialog includes the **Sort by** field. The *Sort by* field displays the name of the column that the Plate List is currently sorted by. To change the content of the *Sort by* field, click on any of the column headers to sort the list by that column. The column header name appears in the *Sort by* field.

#### Plate Status List

The plate status list displays all plates that have been processed. The system database stores all plate data until the plates are deleted. The plate status list contains five columns:

**Plate ID** - displays the ID of the plate.

Assay Code - displays the assay code of the assay that the plate performs.

Status - displays the status of the plate, such as on board, active, and finished.

finished at - displays the time the system completed the results or when they are expected.

last Event - displays a short description of the last event that the plate has been put through.

Click the button at the top of each column to sort the list by that column. The Sort by field in the information area displays the name of the column that the list is currently sorted by.

There are seven navigation buttons on the right of the plate status list that can be used to move through the list. The following navigation buttons are available:

Button	Description
x	Go to the beginning of the plate list.
	Go to the previous page of the plate list.
	Go to the previous entry in the plate list.
Ξ	Select this button to mark multiple entries. Using this button, you can print, cancel or delete all marked entries as a group.
•	Go to the next entry in plate list.
*	Go to the next page of the plate list.
T	Go to the end of the plate list.

#### Buttons

In addition to the common **OK** button, the **Print**, **Processing Steps**, **Delete**, **Detail View** and **Cancel Plate** buttons are also available.





Less than 100000 access rights

This **Print** button is only displayed and accessible if the user has 100000 security access rights. This **Print** button is not displayed if the user does not have 100000 security access rights.

This **Print** button is used to resend selected plate results to the data management software if the initial automated attempt at delivery failed. This failed delivery becomes apparent if the **Plate List** shows a given plate as finished without errors, but this plate does not display with results in the **Results** dialog. A printed copy of the results is generated when the selected plate results are resent using the **Print** button and a confirmation of the data being sent is added to the event log in the status bar (logged with the specific plate identification).

### **Processing Steps**

10:30 10:40 10:50	11:00 11:10 11:20 11:30 11:40 11:50 12:00 12:10	>
R06100422		
		+
OK		

Press the **Processing Steps** button to view a detailed graphical representation of the assays running on the instrument since the last initialization and that are scheduled for processing. The software displays each assay horizontally, with the plate name to the left of the multicolored bar. Each line corresponds to a slot in the plate loading tower. The time scale is represented using the sliding scale located at the top of the scheduler window. The current time is represented by the moving vertical line, moving from left to right horizontally.

The following buttons are also available:

Button	Description
<	Press this button to scroll left through the schedule.
>	Press this button to scroll right through the schedule.
+	Press this button to expand the scale displayed (more detail).
-	Press this button to reduce the scale displayed (less detail).
ОК	Press the <b>OK</b> button to exit the scheduler, back to the <i>Plate List</i> .

The graphical representation uses different colors to symbolize different assay stations. The colors used are:

Color	Description
	(Charcoal gray) The start of plate processing is delayed before it was originally due to be picked up from the plate loading tower.
	Plate barcode scanning
	Transport
	Pipetting steps. The left and right pipettors use the same color.
	Incubation at room temperature
	Incubation at 37°C
	Plate washing step
	Centrifuging or shaking of a plate
	Camera reading

#### Delete

Press the **Delete** button to delete a plate and all of its results from the database. The information cannot be retrieved again.



**Note**: It is not possible to delete plates run as part of a maintenance task from the Plate List.

### **Cancel Plate**

Press the **Cancel Plate** button to cancel a plate from the database. Select the plate(s) to be cancelled and then press the **Cancel Plate** button. A confirmation dialog will appear that must be acknowledged using the **Yes** or **No** button before the plate will be cancelled.



### **Detail View**

Press the **Detail View** button to display the *Detail* dialog. The *Detail* dialog displays more detailed information about the plate selected. The *Detail* dialog has six different tabs displaying different information details:

General Information

Assay Protocol

Plate Events

Sample ID

Raw Results

Flags

## **General Information Tab**

ate id. Ni	08404654	
ssay Code:	QC_CcEe	
perator Name:	Stacy	
late started at:	2009-10-29, 10:24:39	
late finished at:	2009-10-29, 10:38:45	
late Type	Nunc (Nunc)	
	no reaction of the factors	

The *General Information* tab displays basic information about the plate. This information contains the plate ID, the assay code, the operator name, date and times when the plate is scheduled to start and finish, the type of the plate, and the manufacturer of the plate. Note that the times are not the times when the plate actually starts and is finished, but the times when it is scheduled to start and finish.

### Assay Protocol Tab



The *Assay Protocol* tab displays the single steps of the assay in a concise form. The *Assay Protocol* tab serves as a short assay description and displays all of the protocol steps as defined when the assay was performed.

## **Plate Events Tab**

2009-10-29, 10:24:32	"2191	12100302	483582 fo	und in 5-la	ine Bay at	track 3, po	os 3''					
2009-10-29, 10:24:32	"0243	10090460	143572 fo	und in 14-	lane Bay a	t track 14,	pos 6"					1
009-10-29, 10:24:32	''0273	38090470	138613 fo	und in 14-	lane Bay a	t track 14,	pos 9"					
009-10-29, 10:24:33	"1273	38090470	110690 fo	und in 14-	lane Bay a	t track 14,	pos 8"					
009-10-29, 10:24:33	"1283	38090470	114302 fo	und in 14-	lane Bay a	t track 14,	pos 7"					
009-10-29, 10:24:33	''1283	38090470	114302 fo	und in 14-	lane Bay a	t track 14,	pos 7"					
009-10-29, 10:24:33	"9363	37100205	479124 fo	und in 14-	lane Bay a	t track 14,	pos 1"					
009-10-29, 10:24:33	"9451	77100041	476331 fo	und in 14-	lane Bay a	t track 14,	pos 2"					
2009-10-29, 10:24:33	"9540	10100113	700179 fo	und in 14-	lane Bay a	t track 14,	pos 3"					
2009-10-29, 10:24:33	"9640	19110109	014084 fo	und in 14-	lane Bay a	t track 14,	pos 4"					
2009-10-29, 10:24:47	Barco	deReader	- Barcode	NU08404	4654							
2009-10-29, 10:24:59	plate r	noved to F	Reader									
2009-10-29, 10:25:11	Camer	a Reader	results sto	red in file "	D:\Galileo'	images\N	U0840465	54_4ae9a5	5c6_0.bbx			
2009-10-29, 10:25:11	StripE:	dist (Read	er)- Result	s:								
2009-10-29, 10:25:11	66	64	65	65	65	65	64	65	64	64	64	66
2009-10-29, 10:25:11	67	64	64	64	65	65	65	66	66	64	66	67
2009-10-29, 10:25:11	66	64	64	65	65	64	64	65	65	64	64	66
2009-10-29, 10:25:11	65	64	63	64	63	64	64	65	64	64	65	66
2009-10-29, 10:25:11	65	65	65	64	65	64	65	65	64	65	64	67
2009-10-29, 10:25:11	65	66	65	65	66	65	65	65	65	65	65	66
2009-10-29, 10:25:11	64	64	65	64	66	65	66	66	65	65	66	67
2009-10-29, 10:25:11	64	64	65	64	66	65	66	66	65	65	66	67
												F

The *Plate Events* tab displays a log of all events the plate has been through. The event log displays the time each event took place and specific information about the modules where the event took place, for example, *Pipettor: Not enough Liquid LISS*. The plate events are not the planned or expected stations in the assay (as displayed under the *Assay Protocol* tab), but are recorded logs of what actually happened to the plate at a particular time.



**Limitation**: The timestamps for assay plate events that occur simultaneously with the end of Daylight Saving Time (DST) may not accurately reflect the actual time of those events, and it is not possible to predict what these timestamps will display due to their erratic nature. However, in this instance, the assay is successfully completed and the results are unaffected. These erratic timestamps are limited to the end of DST. By running the instrument initialization after the end of DST, the time stamping process is corrected. The timestamps for assay plate events that occur simultaneously with the beginning of DST are accurate and do not exhibit the same erratic features as the timestamps at the end of DST.

## Sample ID Tab



The *Sample ID* tab displays a grid in which each cell represents a plate well; well A1 is in the top lefthand corner. The tab displays the ID of the sample pipetted into each well in the corresponding grid cell. If a well was used for a control or blank measurement, the system also displays it.

On this and other windows that display a grid, clicking on a cell displays more detailed information on that sample or control in the detail window to the right of the **OK** button. Details include well, sample ID, raw result, and any comments. To the right of the **OK** button is a small, complete image of the plate displayed.

### **Raw Results Tab**



The *Raw Results* tab shows a grid in which each cell represents a plate well; well A1 is in the top lefthand corner. The tab displays the raw values from the analysis of the pictures made by the camera reader. The values vary between -2 and +100, where a negative value is an error code. A result marked in red is flagged. You can see the flag when you press on the results well image.

## Flags Tab

	1	2	3	4	5	6	7	8	9	10	11	12
A	2	Reader										
3		Reader										
		Reader										
)		Reader										
		Reader										
		Reader										
		Reader										
1		Reader	Reader	Reader	Reader		Reader	Reader	Reader	Reader	Reader	Reader

The *Flags* tab displays a grid in which each cell represents a plate well; well A1 is in the top left-hand corner. The tab displays all flags assigned to the measured result in the corresponding grid cell. The tab displays a message in the large field at the bottom that explains the specific error message when you click the flag in the cell. This method explains all error flags.

## **Test Results**

#### **About Test Results**

To view test results, click the **Results** button.



The system displays the *Results* screen, with the most recently used screen view (sample or plate).

From the *Results* screen, you can:

View Test Results in Sample View

View Test Results in Plate View

Refer to **Chapter 7 – Test Results** for more information regarding test results.

#### **Before You Begin**



**Limitation**: Validate Method is otherwise known as reflex testing. Reflex testing can only be ordered once per assay in response to results for a given sample number. The instrument WILL NOT order reflex testing again if that sample (with the same barcode) is repeated with the same assay on the instrument at a later phlebotomy date. This may cause incomplete test results to be indicated as final and made available for export.

For example, when a sample is run and assigned the result of A Pending, the reflex testing procedure orders the reflex assay of Weak\_D. If the sample tests positive for Weak D, then the final result of A Weak D will be released. A few days later, another sample with the same barcode is tested on the instrument and the result of A Pending is assigned again. However, on this occasion, the instrument WILL NOT order the reflex assay of Weak\_D, and this second result of A Pending will be the final result. The historic and latest final results are not identical.

This limitation also applies to the reflexive ordering of antibody ID testing for antibody screen positive samples. The consequence is that the second antibody identification test in the history of a given sample number will not be ordered for the second antibody screen positive sample, and so on for subsequent samples.

After the instrument database is archived and the data is removed from the database, the reflex testing procedure will again order a reflex assay on a given sample ID that requires it, because the current sample ID is again new to the instrument software.



**Limitation**: After the reflex testing procedure begins on a sample ID, the value of the result element not being resolved is not updated. For example, if a sample ID is run using the ReflexABO assay and is assigned the result of B NTD, the reflex testing procedure will order a repeat of the ReflexABO assay to resolve result element 2 (the Rh result). If the result from this repeated ReflexABO assay is O Rh Positive, the reflex testing procedure will continue and the current result would be B Rh Positive. The ABO result of O that was determined during the repeat testing is not considered for incorporation into the final result.



**Note**: You must review any test results that have an associated warning flag (<sup>(1)</sup>) through Results. This is to determine if any causes for the warning flag justify any instrument corrective actions. Consultation with Technical Support may be required.

## Maintenance

## **About Maintenance**

The *Maintenance* dialog allows you to view and activate automated maintenance actions. Click the **Maintenance** button.



Refer to **Chapter 10 – Maintaining the NEO Iris** for more information regarding the *Maintenance* dialog.

## Utilities

## **About Utilities**

The *Utilities* dialog allows you to view event logs and statistics, as well as archive and print reports. Click the **Utilities** button.



The system displays the Utilities dialog.

Eve	nt log	Archive Statistics Reports		Utilities		ОК
		Source	Assay Name	Result	File Name	
	1	UD1096524	OCTEST	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
	2	ID02100954	Ab ID	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	Conv
	3	READER	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
	4	G07600103	4 Cell	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
	5	RBC15	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
	6	G07600103	4 Cell	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	<b>D</b> : 1
100	7	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	Print
l ≣	8	PROGRAM2	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
5	9	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
Į Ž	10	UD1095613	QCTEST	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
l Ş.	11	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
Ę.	12	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
E I	13	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
1	14	10-1RVERIFY	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
l ä	15	10-1RPERFORM	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
<u> </u>	16	10-1RPERFORM	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
Į į	17	10-1RPERFORM	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
μū	18	10-1RPERFORM	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
E	19	UD1096555	IS_XMatch	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
8	20	UD1096559	IS_XMatch	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
<u> </u>	21	C2300976	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
	22	C2301058	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
1	23	C2301138	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
E E	24	C2301042	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
E E	25	C2301143	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
l H	26	READ	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
E I	27	RUNGORIGASPCMV	CMV	Successful Conversion	D-1MPA1Data14urora1Lipload1BBX-21	
<u> </u>	뜨					

The Utilities dialog is divided into four tabs:

Event Log

Archive

Statistics

Reports



<u>Note</u>: There may be a delay when entering the Utilities dialog before the screen is populated. This is related to the number of entries in the event log. Wait until the screen fully loads before selecting another tab.

## **Event Log tab**

Eve	nt log	Archive Statistics Reports		Utilities		ОК
		Source	Assay Name	Result	File Name	
	1	UD1096524	QCTEST	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
	2	ID02100954	Ab ID	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	Conv
	3	READER	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
	4	G07600103	4_Cell	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
	5	RBC15	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
	6	G07600103	4_Cell	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	Dwinh
100	7	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	Princ
₽	8	PROGRAM2	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
12	9	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
ΙŦ	10	UD1095613	QCTEST	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
18	11	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
2	12	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
	13	C2301061	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
5	14	10-1RVERIFY	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
8	15	10-1RPERFORM	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
<u> </u>	16	10-1RPERFORM	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
1 2	17	10-1RPERFORM	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
1 2	18	10-1RPERFORM	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
15	19	UD1096555	IS_XMatch	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
8	20	UD1096559	IS_XMatch	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
a di	21	C2300976	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
	22	C2301058	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
Ξ.	23	C2301138	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
1 2	24	C2301042	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
E I	25	C2301143	CMV	Successful Conversion	D:\MPA\Data\Aurora\Upload\BBX-2(	
l đ	26	READ	(Unknown)	No Assay Name Match	D:\MPA\Data\Aurora\Upload\BBX-2(	
d d	27	RUNGORIGASPOMV	CMV	Successful Conversion	D-IMPAIDatalAuroralUnloadiBBX-2	
<u> </u>	<u> </u>					

The Event Log allows you to view the event log, copy it to a disk, or print it.

The *Event Log* contains a list box with a separate vertical tab for each log (**Result and Worklist files**, **Program Error** and **Operator Activity**). The system populates the list box with the contents of the log you select. Initially, no lines are selected. If you click on a line to highlight it (blue) and then click **Copy** or **Print**, only the selected lines are copied or printed.

Column headers are displayed that describe various facets of the event log entries. The column headers are shown below.

Source	Assay Name	Result	
File Name		Created Date	

The window of the dialog displays the event log entries. The order of these entries can be changed if one of the header descriptions is preferred as a sorting mechanism. For example, by pressing the **Assay Name** header, the order of all of the entries in the window will be rearranged according to alphanumeric assay name order.

### Archive tab

Refer to **Chapter 10 – Maintaining the NEO Iris** for information regarding the creation of archives.



**Limitation**: In the event of an instrument computer hard drive crash, all data will be lost from the computer. Any data that was not previously archived will be lost.

#### **Statistics tab**

Event log Archive	e Statistics Reports	OK
Report Name : We Sar We Sar Nur Disl Disl Disl	Il Reaction Results mple Interpretation Results mple Interpretation Result Edits merical Range of Result Edits merical Range of Results tribution of Reaction Values by Date tribution of Reaction Values by Lot Number sts per Assay	Provides a report on the number of tests performed on an Assay within the parameters selected.
Assay :	All Assays	Current Results Only
Operator :	All Operators	1
C Interpreted	All Results Element 1 Element 2 All Results	Postor         1         2         3         4         5         6         7         8         9         10         11         12           A         O
Dates		
Date Search :	All Dates 💌	
A	C Landscape	Clear Test Data Print



**Limitation:** Other than the Test per Assay Report, the outputs from the *Statistics* tab within the Utilities dialog must be verified on site by the end user. Immucor makes no claims as to the accuracy of those outputs.

The *Statistics* tab enables you to select from and customize the available statistical reports. You can customize the reports using the *Filters* options. Only data stored on the NEO Iris computer that remains after any previous archiving procedures can be used as a resource for generating statistical reports.

The **Report Name** list box contains the available statistical reports. The list includes **Tests per Assay** which prints a summary of the number of tests processed per assay. The list is:

- Well Reaction Results
- o Sample Interpretation Results
- o Well Reaction Result Edits
- o Sample Interpretation Result Edits
- o Numerical Range of Results
- o Distribution of Reaction Values by Date
- o Distribution of Reaction Values by Lot Number
- o Tests per Assay

Whichever of the report names is selected from the *Report Name* list, a brief description of that report is displayed in an open area of text to the right of the report list and below the **OK** button. The example shown below is for the *Well Reaction Results* report.

Provides a report on the number and percentage of Well Reaction Results of the given type for the given assay.

The **Assay** drop-down list contains each available assay name, including *All Assays*. You can indicate which assay type to include in the report. *All Assays* is the default.

The **Operator** drop-down list contains each available operator ID, including *All Operators*, which is the default.

The **Results** area has two (2) options, *Well* and *Interpreted*. Each option is only available if an appropriate type of report is selected from the *Report Name* list. For example, if the *Well Reaction Results* report is selected, then both options are made available for use. If the *Sample Interpretation Result Edits* report is selected, then only the *Interpreted* option is made available for use. If the *Tests per Assay* report is selected, then neither option is available for use.

For each option, there is a drop-down list. *Interpreted* has two (2) drop-down lists, *Element 1* and *Element 2* that contain each pre defined result type, including *All Results*. This allows you to specify the types of results to include in the report. *All Results* is the default.

C Well	All Results		-
	Element 1	Element	2
Interpreted	All Results	All Results	•
Dates	All Results NTD		

*Well* has one (1) drop-down list that contains defined well reaction categories, including *All Results*. This allows you to specify the types of well results to include in the report. *All Results* is the default.

- Results		
<ul> <li>Well</li> </ul>	All Results	
	All Results	
	Positive	ľ
s incorproced	Negative	
- Datas	Equivocal	

The **Position** area contains a plate map graphic. The system displays selected wells with a green background and a check mark and non-selected wells with a yellow background and no check mark. Click on a well to toggle it from available to unavailable. Click on a row letter or column number to toggle the entire row or column. Click the *Select All* indicator ( in the upper left of the plate graphic to select the entire plate. These settings indicate which wells to include in the report.

The **Date Search** drop-down list allows you to specify dates for data isolation. The default is *All Dates*. Other options are *Year, Month, Day, Year to Date* and *Date From/To*. Some options have additional date filtering mechanisms that appear in the *Dates* area. The *All Dates* and *Year to Date* options have no additional date filtering mechanisms. If you select the *Year* option, an additional list is made available to allow you to select a specific year.

Dates	2	
Date Search :	Year	<u>•</u>
	2010	

If you select the *Month* option, an additional drop-down list is made available to allow you to select a specific month.

Dates Date Search :	Month	
	December	•
Report Orientatio	January February March April May June July August	

If you select the *Day* option, an additional calendar is made available to allow you to select a specific day.

	2010	/03/0	)4		-	1			
Report Orientation	•	1	Mar	ch 2	:010		F		
A) 2	Sun	Mon	Tue	Wed	Thu	Fri	Sat		
•	28	1	2	3 10	4	12	ь 13		
	14	15	16	17	18	19	20		
	21	22	23	24	25	26	27		
	28	29	30	31	1	2	3		
: VENDORSERVIC	4	5	6	7	8	9	10	000	

If you select the *Date From/To* option, two (2) additional calendars are made available to allow you to set up a date range. If a *From* date is selected from the calendar on the left that occurs after a selected *To* date from the calendar on the right, then both dates are automatically adjusted to reflect a logical date range, which may not be the range you intended to use.

Date Search ;	Date From / 10								
	2010/03/01	2010/03/03 💌							
Report Orientati	on	• March 2010							
	C Portrait	Sun	Mon	Tue	Wed	Jhu	Fri	Sat	
A	• Landscape	28	1	2	3	Ð	5	6	
		7	8	9	10	11	12	13	
		14	15	16	17	18	19	20	
		21	22	23	24	25	26	27	
		28	29	30	31	1	2	З	
	ACE is seen of the worlds seen as he would be	1	-	1.0	-				

**Report Orientation** allows most statistical reports to be printed as either **Portrait** or **Landscape** format by selecting either of the two radio buttons. For some reports, one of the two radio buttons is preselected according to the report name chosen.

The **Current Results Only** checkbox determines if the report is only populated with current results. All report selections, with the exception of **Sample Interpretation Results**, have the checkbox configuration fixed. The **Current Results Only** checkbox for the **Sample Interpretation Results** option can be checked on or off.

The **Clear Test Data** button, located next to the **Print** button, is only visible and available when the **Tests per Assay** report is selected and if you have access to the button within your security setup.



The **Clear Test Data** button is used to delete data because this database is not cleared during the archive process. Clearing the test data should be performed at a regular interval after the data has been analyzed. Refer to **Chapter 10 – Maintaining the NEO Iris** for details of this maintenance task.



The *Tests per Assay Report* prints a summary of the number of tests processed per assay based on the search range chosen. You can narrow the search by assay, operator and/or date range. This report prints a summary of the number of tests processed per assay. The database that contains these details is separate from the NEO Iris result databases and the data remains in the database after archiving. This enables the *Tests per Assay Report* statistics to be printed independent of the scheduled archiving tasks. A sample report is shown below.

Tes	ts per Assay Report				
Printed: 2010/03/02 14:18:20					
Printed By: Vendorservice		<u>12</u>			
Search Parameters					
Assay: All Assays					
Operator: Yolanda					
Date: All Dates					
Assay	Number of Tests	Serial #			
2 Cell	0				
3 Cell	0				
4_Cell	0				
Ab ID	0				
ABO AB 2	0				
ABORH	0				
ABORH 2	0				
ABORH_AB	0				
Ag_C RH2	0				
Ag_c RH4	0				
Ag_CcE	0				
Ag_CcEe	0				
Ag_E RH3	0				
Ag_e RH5	0				
### **Printing Statistics**

To print statistics:

Step	Action
1	After you select the report and filters, and then press the <b>Print</b> button, the system displays a preview window with a progress bar showing each assay, one after another.
	Well Reaction Results - Statistics Duppet         Augent 2         Augent 2         Cancel
2	When the system completes the statistical report, you can view the report before printing.
	Sample Interpretation Results - Statistics Durput         Image: Interpretation Results - Statistics Result           Sample Interpretation Results - Statistics Result         Result - Statistics Result           Image: Results - Statistics Result         Result - Statistics Result           Image: Results - Statistics Result         Result - Statistics Result           Image: Results - Statistics Result         Results - Statistics Results           Image: Results - Statistics Results         Results - Statistics Results           Image: Results - Statistics Results         Results - Statistics Results           Image: Results - Statistics Results         Results - Statistics Results           Image: Results - Statistics Results         Results - Statistics Results           Image: Results - Statistics Results         Results - Statistics Results           Image: Results - Statistics Results         Results - Statistics Results           Image: Results - Statistics Results         Results - Statistics Results           Image: Results - Results - Results         Results - Results           Image: Results - Results - Results - Results         Results - Results - Results           Image: Results - Results
	With the slider bar in the bottom left corner you can enlarge or minimize the preview. Click <b>Print</b> to display the system print dialog and print the report.

#### **Reports tab**

The **Reports** tab is an area to configure the upload pathways of instrument report formats. This area is only available with the correct user access rights. Refer to **Chapter 8 – NEO Iris Reports** for details of report formats.

### **Instrument Settings**

#### **About Instrument Settings**

The *Instrument Settings* dialog allows you to program the system settings in order to configure the instrument to the individual needs of the laboratory. Press the **Settings** button.



The system displays the Instrument Settings dialog.

	5
Setup           Assay Definition         Profile Definition         Setup Barcodes         Setup Device         Setup Plates         Users         Maintenance         Service         LIS         Security	ОК
Assay List	
B         Reflext/WD           B         ResVol           B         ResVol           B         Rev_AB0           B         Reform           R         Perform	
B         ServiceF1           B         Syph           B         Syph, RPT           B         SyringeEx           B         Weak_D           B         Weak_D_F	
Edit New Save	



**Note**: The Settings button is not available for access when the instrument is actively processing samples and plates or is initializing, regardless of your security access level. This prevents the changing of settings during instrument activity.

Access to different tabs is controlled by user access rights and some tabs are limited to Immucor personnel only. Operator accessible functions are described below.

From this dialog, you can program user names, passwords, and access rights. For complete instructions, refer to **Assigning Passwords and User Access Rights** in Chapter 4 - Security.

#### **Barcode Settings Configuration**

Under the Setup Barcodes tab of the Setup dialog, you can configure barcode features by checking or clearing check boxes or radio buttons for both Donor and Sample barcode information. By checking a box or selecting a radio button, you make that particular feature active.

	Setup
Assay Definition   Profile Definition   Setup Barcodes	Setup Device Setup Plates Users Maintenance Service LIS Security
Donor Sample Mask: Prefix: Codabar Conversion G ABC Custom ISBT 128 Conversion	ID Mask Directions: Use wildcard character(s) to represent expected ID mask. ? Ignore character # Accept character (numeric only) & Accept character (alpha or numeric) * Accept tremaining characters (alpha or numeric) Use a numeric range to select the wanted characters from an ID. Examples N9835982 using the mask ??* would result in 835982 as the ID. 001FG37825d2 using the mask ???########? would result in 37825 as the ID. 001FG37825d2 using the mask ???########? would result in 7637825 as the ID. A2140119KS using the mask 2-8 would result in 2140119 as the ID.
Assay Definition     Profile Definition     Setup Barcodes     Setup D       Donor     Sample       Mask:	Assay Definition Profile Definition Setup Barcodes Setup D Donor Sample Mask: Codabar Conversion: C ABC C Ogston
Donor	Sample

**Mask:** Is used to define how a barcode is manipulated and is governed by how characters from the legend are combined. Masking is an alternative flexible method to manipulate the barcodes.

? means to Ignore a barcode character.

# (Digit) means to use a barcode character if it is a digit.

& means to use Any character.

\* means to use All remaining characters.

2–9 is an example of a numeric range that defines which characters are to be used. This number range can be customized.

The following ISBT 128 barcode data is an example of masking. If the mask was set at ?2-14, then the ISBT 128 barcode of =G07359931656200 would be masked to G073599316562. Additional examples are listed in the explanatory text displayed on the Setup Barcodes tab.

You can add a prefix to the barcode data by entering the desired prefix into the **Prefix:** field. The maximum number of characters that can be entered into the prefix field is five. **Prefix** is only available on the **Donor** tab.



**Note:** Only one prefix per rack can be added at any one time.

**Warning:** You must be careful to make sure that the samples in one given rack are of a single prefix type. It is your responsibility to check this data. If you have mixed samples of different prefixes, all samples will have the same prefix added to that which you entered. Thereby, identification errors could occur if the operator fails in this tube data check.

The Codabar Conversion feature activates the alphanumeric conversion for Codabar barcodes. In the case of ABC conversion, the first two digits (of seven) of the barcode are replaced with designated alphas using a software embedded truth table. The Custom conversion is based on a site configurable software embedded truth table.

The ISBT 128 Conversion feature activates the stripping of the equal sign (=) at the beginning of the sixteen-character ISBT 128 blood unit identification barcodes and also the last two checksum digits, thereby leaving characters numbered between two and fourteen.

Even though settings can be configured, you can temporarily edit these settings on a rack-by-rack basis using the Prefix + Barcode Settings button on the *14-lane Bay* dialog.



When exiting and subsequently reentering the *14-lane Bay* dialog, the pre-configured settings will again be the default. The pre-configured settings cannot be changed by temporarily changing the settings on a rack-by-rack basis.

If a donor prefix is configured, it will be displayed as replacement text on the *14-lane Bay* dialog *Prefix + Barcode Settings* button. For example, if a donor prefix of *12345* is used, then it will display, as shown below, as replacement text on the button.

r Prefix: 12345	
	INFIENX 12343

Refer to Chapter 6 – Instrument Testing Operation for more information regarding Barcode Settings.

You can close Setup by pressing the **OK** button at the top of the dialog to save the settings.

# **Creating Assay Profiles**

You can use the *Profile Definition* tab to create assay profiles:

Assay Definition P	rofile Definition	Setup Barcodes Setup De	vice Setup P
Profile		Assay           2_Cell           3_Cell           4_Cell           ABORH_2           AQ_CRH2           AQ_CRH2           AQ_CRH4           AQ_CRH4           AQ_CRH5           AQ_CR15           AQ_CR15           AQ_CR15           AQ_CR15           AQ_Q00           Flow3	*
Default Profile: New Rename	Delete Default	Add	Select All

Step	Action
1	Press the <b>Settings</b> button from the Main Menu Bar.
2	Select the Profile Definition tab.
3	Press the <b>New</b> button from the Profile box.
4	Type a profile name into the open text box with the flashing cursor, for the new profile to be created.Image: Mote: Greater than twelve (12) characters can be typed into the open text box, however only a maximum of twelve (12) characters can be displayed on the assay profile button. If greater than twelve (12) characters are typed into the open text box, then the excess characters will not be visible on the profile button.
5	Select the assay to be included in the profile from the Assay box list.
6	Press the <b>Add</b> button from the Assay box. <b>Note</b> : The <b>Select All</b> and <b>Unselect All</b> buttons can be used to collect all assays into a profile or remove all assays from a profile respectively.

7	Repeat steps 5 and 6 for any additional assays to be added to the profile.
	Note:If the newly created assay profile is required to be the default assay profile assigned to all samples tested on the instrument, then you must press the Default button after step 7 of this procedure. The system displays the assay profile name next to the Default Profile: text. You can remove the default assay profile status by pressing the Default button a second time, prior to step 8 of this procedure. You can only select one (1) assay profile as the default assay profile at any time.
8	Press the <b>OK</b> button on the <i>Setup</i> dialog when you are finished.
9	You will be asked to enter your login password in an Operation confirmation dialog. After your password entry, press the <b>OK</b> button of the Operation confirmation dialog.

The **Rename** button can be used to rename an existing profile name.

# Help

#### **About Help**

The *Help* window displays the electronic NEO Iris Operator Manual on the NEO Iris PC monitor as well as providing access to the remote support functions and listing module version information.

#### Accessing Help

Press the Help button on the Main Menu Bar to display the Help window.



You must insert the current CD version of the electronic operator manual into the CD drive of the instrument PC tower to provide the files for the monitor display of the manual.



**Note**: Ensure that you have loaded the current CD version of the electronic NEO Iris Operator Manual into the CD drive of the instrument PC tower before pressing the Help button. If the CD drive is empty, or contains a wrong or damaged CD, the system displays an error message indicating that the CD of the manual must be loaded into the CD drive.

mpa		×
4	Cannot find the "e:\NEO_Help.chm" file. Check to see that the file exists on your disk. If it doesn't, you need to reinstall it.	
	OK	

The system displays the *Help* window.

	Context Help (12023)
Extended Help	
Version Info	
Remote Support	
ОК	

#### **Extended Help**

When you press the **Extended Help** button (only visible when the NEO Iris Operator Manual CD is loaded), the software displays the *NEO Iris Operator Manual* window at the **About this Manual** section.

Import of the formation of the for	About this Manual In this section In this section In this section In the section of the section
Concerning Concerning Technical Disa Signature Concerning Concerning Technical Concerning Concerning Manual Concerning Concerning Manual	

Press the hyperlinked text of a topic to access the selected topic. The display of the help content list is dependent on whether it was displayed when the file was last accessed on the PC. If the list was hidden on the most recent access, when the help display was closed at that time, then it will not display, or the list will be shown if it was displayed when on-line help was closed at that time. You can then press the *Hide* or *Show* buttons to toggle the display of the help contents. Pressing on any of the content topic links will display the selected topic in the help window.

Click the **X** button at the top right of the *NEO Iris Operator Manual* window to exit the Extended Help page and return to the Help window.

#### **Version Information**

Click the **Version Info** button. The *Version Information* dialog is displayed.

Version Information (Instrument 5/N 5030090010)		
	⊡- NEO Iris: 1.6.11.1 (2016-04-05)	<u> </u>
	COP	
	<ul> <li>DataReduction (bxdr.dll): 1.6.11.0 (2015-11-10)</li> </ul>	
	IOModuleAdapter (ioma.dll): 1.6.11.0 (2015-11-10)	
	BarcodeReader (bcrd.dll): 1.6.11.0 (2015-11-10)	
	🚊 · Reader (rdis.dll): 1.6.11.0 (2015-11-10)	
	IStar Dynamic Link Library: 1.8.4.1	
	IStar Interface: 1.8.1.0	
	🚊 · Centrifuge (cefu.dll): 1.6.11.0 (2015-11-10)	
	Firmware: 04.0006 5032 2010-11-10	
	🚊 Washer1 (wash.dll): 1.6.11.0 (2015-11-10)	
	Firmware: 04.0019 2274 2013-02-08	
	🖻 Incubator1 (incu.dll): 1.6.11.0 (2015-11-10)	
	Firmware: 03.0003 2255 2010-05-03	
	Pipettor (pip4.dll): 1.6.11.0 (2015-11-10)	
	🚍 Transport (pdrv.dll): 1.6.11.0 (2015-11-10)	
	È Firmware: 01.0018 2251 XDV 2012-04-13	
	z-Axis: 01.0011 2251 ZDV 2015-07-24	
	y-Axis: 01.0009 2251 YDV 2015-07-24	
	🖻 Sample Loading Bay (pald.dll): 1.6.11.0 (2015-11-10)	
	⊟ Firmware: 04.0018 4292 LBY	
	- 14-lane Bay Scanner: BCL 148.1 V 50.15 23.02.07	
	5-lane Bay Scanner: BCL 21 V 03.58 05.02.2002	
OK	🚍 · Tower (ptwr.dll): 1.6.11.0 (2015-11-10)	-
UN		
	<u>, , , , , , , , , , , , , , , , , , , </u>	

The Version Information dialog displays a list of the modules and software used on the instrument.

The dialog displays the information in a tree structure with modules at the roots and their dependent sub-modules as the branches. Press a plus sign (+) to expand the list and view more information under a topic. Press a minus sign (–) to collapse the list and hide information. Press the **OK** button to exit the *Version Information* dialog.

#### Access to Remote Support

Access to remote technical support using *blud\_direct* is provided through the use of the **Remote Support** button in the help window:

Step	Action
1	Press the <b>Help</b> button without the help CD being loaded into the PC disk drive.
	<b>Result</b> : the following error message will be displayed:
	Cannot find the "e:\NEO_Help.chm" file. Check to see that the file exists on your disk. If it doesn't, you need to reinstall it.
	ОК
	Note: The error message will not be displayed if the help CD is loaded in the PC disk drive. In that instance, skip step 2 of this procedure and proceed to step 3.
2	Press the <b>OK</b> button of the error message.
	OK
	Result: The error message disappears and the help window is displayed.
3	Press the <b>Remote Support</b> button.
	Remote Support
	Result: The blud_direct Enter session: page is displayed.
	Enter session:
	Send

4	You must then enter the correct session ID (provided by Immucor Technical Support over the phone) into the Enter session: field and then press the Send button within fifteen (15) minutes of assignment, otherwise it will expire and a new session ID will have to be assigned to you. The session ID is not case sensitive. When the session ID is confirmed, you will be presented with a message to accept contact by Immucor Technical Support (operator).
5	Click the Accept button to continue with the support session.
6	The chat window of the blud_direct application will then display. You can then interact as necessary through the chat window which includes responding to any requests by Technical Support for remote control of your instrument computer. Remote access to your instrument computer is only provided when you accept a request to do so while in the current chat session.
	Note: In addition to the chat window as a newly displayed window, the Help page will update with a button graphic. Do not press this button because doing so will create a new chat session.

#### blud\_direct options

The table below describes the blud\_direct options.

Option	Description
·Ť· TĐ TĐ	The font of the text displayed within the chat can be increased, decreased, or reset by clicking the appropriate buttons shown.

💽 Receive Chat Transcript	To receive a transcript of the details regarding your support
	session, click on the <b>Receive Chat Transcript</b> button.
	You will be prompted to enter where the transcript should be sent.
	Enter the e-mail address to which you would like to have the
	transcript sent and then click the <b>Accept</b> button.
	C Receive Chat Transcript
	Enter your e-mail address and we will send you a complete transcript of your conversation with the operator
	Your e-mail:
	Accept Cancel
	Note: You may click on Receive Chat Transcript at any time during the chat session, but must do so before closing the chat window (i.e. before ending the session).
	It is important that a valid e-mail account is entered into the prompt. Please remember to check spelling, etc. of the address before clicking the <b>Accept</b> button.
Display chat window again	Should the chat window fall behind the instrument software, close the <i>Help</i> page and then reopen it. Press the <b>Display chat window</b> <b>again</b> link. This will allow the chat window to display again in front of the instrument software.
	Do not minimize the chat window. This will prevent the window from being able to be displayed again using the <b>Display chat</b> window again link.
	To end a remote control session without ending the chat session, you can press the <b>Terminate</b> button. To end the chat session, you can press the close icon in the upper right corner of the chat window. In closing the chat session, the remote control session will also be ended.

## Shutdown

#### **About Shutdown**

The *Shutdown* dialog guides you through the procedure that you must perform before the system is shutdown. Press the **Shutdown** button.



The system displays the *Shutdown* dialog.

Refer to Chapter 9 – System Shutdown for information regarding system shutdown.

# Machine Monitor

### About the Machine Monitor

#### Parts of the Machine Monitor

The Machine Monitor is a graphical view of the instrument deck that displays the status of each of the system modules and provides access to the module dialogs. The instrument displays this screen when the system is in standby mode.



The table below describes the sections of the Machine Monitor.

Area	Name	Description
A	Washer	By pressing the button, the <i>Wash Buffers</i> dialog is displayed. This button also provides an overview of the system liquid on the instrument.
В	Camera Reader	By pressing the button, the most recent successfully captured image taken by the camera is displayed.
С	Incubator	By pressing the button, the <i>Incubator</i> dialog is displayed. This button also provides an overview of the incubator status.

Т

Г

D	Plate Tower	By pressing the button, the <i>Plate Loading Tower</i> dialog is displayed. This button also provides an overview of the plates in the loading tower.
E	14-lane Bay	By pressing the button, the <i>14-lane Bay</i> dialog is displayed. This button also provides an overview of the racks loaded in the 14-lane bay and the status of the module.
F	5-lane Bay	By pressing the button, the <i>5-lane Bay</i> dialog is displayed. This button also provides an overview of the racks loaded in the 5-lane bay and the status of the module.

### Incubator

#### **Incubator Dialog**

The *Incubator* dialog provides an overview of the status of the incubator. Click the **Incubator** button in the Machine Monitor to display the *Incubator* dialog.

Incubator			
	state	set-point [°C]	actual value [°C]
15	free	39	39
14	free	39	39
13	free	39	39
12	free	39	39
11	free	39	39
10	free	39	39
9	free	39	39
8	free	30	35
7	Isolation	20	26
6	free	20	21
5	free	20	23
4	free	20	24
3	free	20	25 Environment
Done 2	free	20	25
1	Plate in Slot	20	25

The Incubator dialog displays the Incubator Status List and the Environment area.

**Incubator Status List** - The incubator status list displays the status of the individual incubator positions. The State column (15–1) represents the positions in the incubator. If a plate is present in a position, it is indicated here. The Set Point column displays the temperature set for this incubator position, such as 39°C. The Actual temperature column displays the actual temperature of this incubator position. The acceptable range of a 39°C incubator position is 38°C to 40°C inclusive. Assay results are invalidated if the temperature falls outside of this range during the assay incubation.

**Environment Area** - The Environment area in the lower right-hand corner of the *Incubator* dialog displays the actual ambient temperature measured inside the instrument.



**Note:** The thermodynamics of the incubator bays require an actual temperature of 39°C to yield a test well temperature of 37°C.

### Wash Buffers

#### Wash Buffers Dialog

The *Wash Buffers* dialog provides an overview of the system liquid and the waste level status. Click the **Washer** button in the Machine Monitor to display this dialog. You can also access the *Wash Buffers* dialog from the *Resource Overview* window.

Wash Buffers	
Fill Level: In use Buffer ID: PBS Information: Assay: 2_Cell Buffer: PBS	
Initial Prime Remove	
Done Vaste	✓

The Fill Level area displays a green check mark ( $\sqrt{}$ ) if the liquid level is above the threshold level, or a red exclamation mark (!) if the liquid level is below the threshold level. The threshold is set in such a way that the remaining volume is sufficient for the completion of the plate currently in the washer.

If you use the **Initial Prime** button at the bottom of the wash buffer container area to request an extra prime for the system liquid, there is a delay of several seconds between pressing the **Initial Prime** button and pump activation.

Refer to **Filling the System Liquid** in **Chapter 10 – Maintaining the NEO Iris** for more information regarding the Wash Buffers dialog.

### Reader

#### **Reader Dialog**

Press the **Camera Reader** button to display the *Reader* dialog.



The *Camera Reader* dialog displays the most recent successfully captured image taken by the camera. The name of the **Assay** and **ID** associated with the plate are also displayed. Press the **Done** button of the *Reader* dialog to return to the *Machine Monitor*.





**Note:** The above image displays the Reader Verification Plate tool.

### **Plate Loading**

#### **Plate Loading Dialog**

The *Plate Loading Tower* dialog provides an overview of the plates in the loading tower and prompts the operator when loading plates. Press the **Plate Tower** button in the Machine Monitor to display this dialog. You can also access the *Plate Loading Tower* dialog from the *Resource Overview* window.

Plate Loading Tower	-	
Scan new plates Cancel Schedule	No Status         Plate ID         Assay Code           15	Strip Selection       Scan Plate       Assay Selection       Expiry Date       Processi • •         Assay Selection       3_Cell       4_Cell       -         Ab_ID       ABORH       ABORH_2       -         Ag_C RH2       Ag_c RH4       Ag_CCEe       -         Ag_E RH3       Ag_e RH5       Ag_Kell       -
Done		



**Note**: Use the left and right arrow buttons at the bottom of the *Assay Selection* tab to scroll through the available assays.

The *Plate Loading Tower* dialog is divided into two areas, the loading tower diagram on the left and the tabs on the right. There are several buttons in the *Plate Loading Tower* dialog that have functions described in the table below.

Button	Description
4 1	Press these buttons to switch the display of tabs in the left or right direction.
Scan new plates	Press the <b>Scan new plates</b> button to initiate the barcode scanning of new plates.

Button	Description	
Cancel Schedule	The <b>Cancel Schedule</b> button on the main <i>Plate Loading Tower</i> dialog is used to cancel all plates that have not begun processing. A <i>Cancel Confirmation</i> will appear that should be acknowledged using the <b>Yes</b> or <b>No</b> button before all plates that have not started processing will be cancelled. Only plates that have not started processing can be cancelled. Such cancelled plates retain the assays scheduled and will reappear in the Resource Overview window. If these plates are cancelled, any unused reagent will be electronically reassigned back to its vial. Such reagent reassignment is also activated if these plates are aborted during error recovery and troubleshooting.	
Done	Press the <b>Done</b> button to exit the <i>Plate loading Tower</i> dialog.	

#### Loading Tower Diagram

The loading tower diagram displays the positions in the loading tower (15–1). For each position, the dialog displays the status, the plate ID (barcode), and the codes of the assays that can be performed with this plate. If the plate is already scheduled, the dialog displays only the assay that actually runs with this plate.

#### Tabs

Five tabs are available in the *Plate Loading Tower* dialog: Strip Selection Scan Plate Assay Selection Expiry Date Processing Steps

#### Automatic Barcode Plate Identification

Place the plates into the loading tower with their barcodes facing to the left. For more information, refer to **Loading Plates** in **Chapter 6 – Instrument Testing Operation**. Close the plate tower door after the LED for that slot has changed to orange. Barcodes will be scanned automatically. If the automatic scanning process is not initiated because the plate tower dialog is not closed, press the **Scan new plates** button instead. The transport mechanism removes each new plate and presents its barcode to the laser scanner. The system updates the loading tower graphic to show the barcode of each plate and the assays that it may be used for.

#### Manual Plate Identification – Scan Plate Tab

To load a plate with an unreadable barcode ID, you can type the barcode manually.

Step	Action	
1	Click the <i>Scan Plate</i> tab.	
2	The scan plate screen appears.	
	1. Scan plate ID:	
	2. Put plate in Plate Loading Tower.	
	Type the barcode into the field either manually or using the hand-held barcode scanner. You must enter the plate lot number followed by the serial number if manual entry is used.	
3	Load the plate into a free slot in the loading tower (indicated by a continuous green LED).	
	The instrument allocates the barcode ID to the newly inserted plate. The system uses the plate as if the integrated barcode scanner had read the barcode.	

#### Manual Plate Identification – Assay Selection Tab

You can also use the Assay Selection tab to manually type a barcode and assign an assay to a plate.

Step	Action
1	Click the Assay Selection tab to enter the plate details manually.
2	Insert the plate into an empty slot in the loading tower. The system updates the loading tower graphic.

Step	Action
3	Type the plate ID in the respective position in the loading tower diagram (double-blind entry required). The corresponding assay assignment will automatically be populated if entering a plate barcode; if manually entering an id to run a maintenance task, click the assay button that should be run with that plate.

#### **De-selection of Absent Strips – Strip Selection Tab**

The instrument supports de-selection of strips that are absent from the plate. This feature only works when an assay is already assigned to a plate. By default, the instrument assumes that when you insert a new plate, all strips are present on the plate. If this is not the case, you can deselect strips and make them unavailable for the instrument.



**Limitation**: The camera reader will verify the presence of clean and correctly positioned strips on a plate before starting an assay. If the reader detects the presence of strips that should have been deselected, such as wrong position of used strips, then the plate is aborted.

Select the *Strip Selection* tab and then select a plate.

- Plate Loading Towe	Ma Status - Plata ID	Assau Cada	Strip Selection Scan Plate Assay Selection Expiry Date Processi
Scan new plates			Strip Selection
Cancel Schedule	14		
	13 O	ABORH ABORH_2	
	11		
	10		
	8	-	
	7		
	6 <b>()</b> 5 <b>()</b>	_	
	4		Plate ID: UA4310947
	3		Assay(s): ABORH ABORH_2 FWD_ABORH mAgScr37 Use before: 2030-12-31 On Board Time: 0 days: 0 hours: 4 minutes
Done	1		First loading: 2009-09-10, 11:59:56 QC status: due

By pressing on a vertical strip icon, you can toggle the strip on and off.

If the instrument has used up some but not all strips in an assay run, and the system brings the plate back to the loading tower, you can start the plate again. In this case, however, it is not possible to reuse already used strips for a second time.

#### Assigning Expiration Dates to Plates – Expiry Date Tab

Manual expiry date entry is only required for *Scheme 1* plate frame barcodes and not for *Scheme 2*. If the *Scheme 2* lot number is unknown to the software, such as when the first plate of a new lot is loaded and scanned, the *Expiry Date* tab fields will automatically populate with date information when the plate frame barcode is scanned by the instrument.

When a *Scheme 2* plate barcode is assigned to a loading tower slot, the software disables editing of the plate *Expiry Date* tab.

Refer to **Chapter 1 – Introduction to the NEO Iris** for more information regarding *Scheme 1* and *Scheme 2* plate frame barcodes.

You can assign an expiration date to *Scheme 1* plate frame barcodes. When the barcode scanner scans the plate frame, the barcode is interpreted. If the lot number is unknown to the software, such as when the first plate of a new lot is loaded and scanned, the system displays a message alerting you to enter the expiration date for this plate.

To enter the *Scheme 1* expiration date:

Step	Action
1	Select the <i>Expiry Date</i> tab and then select the plate with a left-trackball click or via the touch screen.
	Plate Loading Tower       No       Status       Plate ID       Assay Code         Scandle       15       Image: Strip Selection       Scan Plate       Assay Selection       Expiry Date       Processi • •         Cancel       12       Ud4310947       ABORH ABORH_2       Plate Loi:       U/A43         10       Image: Selection       Selection       Selection       Selection       Selection         9       Image: Selection       Selection       Selection       Selection       Selection         9       Image: Selection       Selection       Selection       Selection       Selection         9       Image: Selection       Selection       Selection       Selection       Selection         0       Image: Selection       Selection       Selection       Selection       Selection         0       Image: Selection       Selection       Selection       Selection       Selection         0       Image: Selection       Image: Selection       Selection       Selection       Selection         0       Image: Selection       Image: Selection       Image: Selection       Selection       Selection         0       Image: Selection       Image: Selection       Selection       Selection       Sele
	The system displays a new dialog where you can enter the expiration date for this plate lot.

Step	Action
2	Press the <b>Set Date</b> button. The system assigns this date as the expiration date to all plates of
	that lot on the system.
	Stip Selection       Scan Plate       Assay Selection       Expity Date       Processi • •         Use unit       Plate Lot:       12         year       month       day         2030       •       12       •         Set Date       Remove Date       Remove Date
	<b><u>Limitation</u></b> : If an incorrect, but valid, plate expiration date is entered and an

assay using this plate lot is run, then the incorrect expiration date is entered and an assay using this plate lot is run, then the incorrect expiration date is stored for that lot in the software. If the incorrect date is then corrected in the software following this assay run, and this correct date is further into the future than the original incorrect date, the data management system cannot update the incorrect stored date. Consequently, if an assay is run using this plate lot again, at a date of testing between the original incorrect date and the corrected date, then these results will be displayed in the *Plate View* window with an associated warning message of **Plate failed criteria**, **please check Interpretation Failures**. The validity of these results is not compromised by this chain of events that generated this warning message, and these results can still be exported or printed.



/ • \

**Warning**: If an incorrect plate expiration date is entered, which is beyond the actual expiration date of the plate, and an assay using this plate lot is run at a date after the actual expiration date, but before the erroneously entered software expiration date, then the validity of the results is compromised, even though the results will appear legitimate. It is the operator's responsibility to enter the correct plate expiration date into the software.

If you have access to change plate expiration dates and are logged into the software, you can clear the Set Date by pressing the **Remove Date** button. You can then set a different expiration date.

You cannot use a plate without an expiration date. If the expiration date on a plate is exceeded, you cannot use the plate.

#### Viewing Processing Steps – Processing Steps Tab

Select the *Processing Steps* tab to view a graphic of the current schedule.

The scheduler axis line moves in real time from left to right to demonstrate time passing as the instrument operates.



<u>Note</u>: For more information on scheduling, refer to **Chapter 1 – Introduction to the NEO** Iris.

Scan	No Status Plate ID	Assay Code	Scan	Plate A	ssay Sele	ection	Expiry Dat	e Proce	ssing Steps	4
new plates	15	_	×	13:00	14:00	111	15:00	16:00	17:00	1
Cancel	14 🔘	-	×							
Schedule	13 🔵		×							
	12 😝 UA4310947	ABORH ABORH_2	×							
	11 🔘		×							
	10 🔘		×							
	9 🔘		×							
	8 💽		X							
	7 🔘		×							
	6 🕒		×							
	5		×							
	4		X							
	3		×							
1	2		×							100
Done	1 🔵		X							J.

In addition to the buttons on the *Plate Loading Tower* dialog, the following navigation buttons are also available:

Button	Description
<	Press this button to scroll left through the schedule.
>	Press this button to scroll right through the schedule.
+	Press this button to expand the scheduler time scale displayed (more detail).
	Press this button to reduce the scheduler time scale displayed (less detail).

Button	Description										
×	Each plate loading tower position on the scheduler has a cancel button, so that any individual plate being processed can be cancelled. Select the <b>X</b> that corresponds to the plate that should be cancelled/aborted. The specific <b>X</b> button corresponds to the chosen plate that is horizontally level with the button. A <i>Cancel Confirmation</i> will appear that should be acknowledged using the <b>Yes</b> or <b>No</b> button before the selecte plate will be cancelled. Such a cancelled plate retains the assays scheduled and will reappear in the Resource Overview window. If the plate is cancelled before it has begun processing, any unused reagent will be electronically reassigned back to its vi Such reagent reassignment is also activated if a plate is aborted during error recover and troubleshooting.										
	Cancel Confirmation  Do you really want to cancel processing of plate 'R06100422' with assay 'QC3_Cell'?  Yes No										
	<b>Note</b> : Cancelled plates are not lost to the system. Such cancelled plates are subsequently available after the STAT plate is started processing.										

The scheduler graphical representation uses different colors to symbolize different assay stations. The colors used are:

Color	Description
	(Charcoal gray) The start of plate processing is delayed before it was originally due to be picked up from the plate loading tower.
	Plate barcode scanning
	Transport
	Pipetting steps. The left and right pipettors use the same color.
	Incubation at room temperature
	Incubation at 37°C
	Plate washing step
	Centrifuging or shaking of a plate
	Camera reading

# 14-lane Bay

#### Purpose

The **14-lane Bay** dialog provides an overview of the racks loaded on the 14-lane loading bay. In this dialog, you can view the status of the racks, the samples or reagents in each rack, and the assays requested for the samples. You can also use this dialog to enter a sample ID for a position in a rack and to manually request assays.

1 2 3	3 4 5 6 7 8 9 10 11 12 13 14 <sup>Red</sup>	Rack: Prefix +	Barcode Settings
			4_Cel
			ABORH
			Ag_C RH2
5 5 5	**********		-
			Profile

The system displays the 14 lanes of the loading bay. If a rack is in a lane, a rack is displayed in the lane on the icon. The system displays a green check mark ( $\checkmark$ ) at the bottom of the rack on the Machine Monitor if the rack type is recognized. The system displays a red question mark (?) at the bottom of the rack on the Machine Monitor if the rack type is not recognized. The five-lane reagent loading bay exhibits the same features.

#### Parts of the 14-lane Bay

Press the **14-lane Bay** button in the Machine Monitor to display this dialog. You can also access the *14-lane Bay* dialog from the *Resource Overview* window.

Ra	ck Are	ea —		1.92	19. 11	•	14 10	14 14	14 - 14h	8 10	14 14	14. IV.		Test Selection
1	2	3	4	5	6.	<b>A</b>	8	9	10	11	12	13	14	Rack: Prefix + Barcode Settings 2_Cell 3_Cell
2	0	0	0	0	0	0	0	0	0	0	0	0	0	ABORH ABORH ABORH ABORH2 ABORH2
CL	10	0	0	0	CL	B	CL	10	0	CL	5	CL CL	CL	
											[			Profile
	)one					Reca	11		All					<u><u> </u></u>

- A: Lane button
- **B**: Lane
- C: Assay button
- D: Profile area

The *14-lane Bay* dialog is divided into three areas: the Rack Area, the Identifiers area, and the Test Selection area.

#### **Rack Area**

The *Rack Area* displays a graphical, aerial view of the sample loading bay and the **Recall** and **All** buttons. The 14-lane bay graphic displays each of the 14 lanes of the bay with a corresponding button above. If a rack is present in a lane, the dialog displays a graphic of a rack in the lane. Click a rack to select it. A dark gray bar highlights the currently selected rack.

Use a lane selection button (1–14) to focus the barcode reader on the desired lane. The colored bar in the button displays the status of the lane in accordance with the consistent color codes. The color displayed here is the same as that of the lane LED on the instrument.

If a sample rack is loaded and one or more barcodes on that rack is/are unreadable, you can correct the unread barcode. Click the **Recall** button to manually enter unreadable barcodes. Refer to **Chapter 6** – **Instrument Testing Operation** for more information regarding how to rescan unread barcodes.

#### **Identifiers** Area

The *Identifiers* area displays the details of the selected rack of the 14-lane bay graphic. The system displays the rack type in the Rack type field at the top of the Identifiers area, and sample information in the list below.

The first column of the list displays the position in the rack. The second column contains tags that are used to select the type of material in the container. The following options are available:



The Sample tag indicates that sample material is loaded. The Donor and Control tags indicate that donor or control material is in this position. Reag1, Reag2, and Reag4 indicate that reagents are loaded in this position. The number states the number of probes of the 4-probe arm that can access the vial simultaneously. The tags are not strict. The material in a vial labeled Sample can also serve as a donor. Controls can also be labeled as reagents. The tags are set by default when you insert a rack. To change them manually, you must activate the ID field, and then press the button to change the tag.

If the rack contains samples with barcodes, the system displays these in the text box of the third column. If a sample does not have a barcode, or the instrument is unable to read the barcode, you can enter the ID manually. Refer to **Chapter 6 – Instrument Testing Operation** for information regarding the input of barcodes that are not initially scanned by the instrument.

The **Prefix + Barcode Settings** button can be used if barcodes require any data conversion or manipulations. Refer to **Chapter 6 – Instrument testing Operation** for information regarding the use of this button.

#### **Test Selection Area**

The *Test Selection* area contains two sets of buttons. At the top are buttons for all assays available on the system. At the bottom are buttons for all of the profiles (predefined combinations of assays) defined on the system. Use these buttons to manually assign an assay or profile to a sample.

To assign an assay or profile, click a sample ID in the *Identifiers* area and then click an assay or profile button. If the system does not display an assay or profile button in the Test Selection area, use the **Previous** and **Next** buttons to scroll through the pages of the assay list.

The **STAT** check box can be used if samples must be tested urgently. Refer to **Chapter 6 – Instrument Testing Operation** for information regarding the use of the **STAT** feature. If a reagent rack is loaded or selected for display, the software operates as described in 5-lane Bay below.

# 5-lane Bay

#### Purpose

The **5-lane Bay** dialog provides an overview of the racks loaded in the 5-lane bay. This bay is designed for the use of reagents only. In this dialog, you can view the status of the racks and the reagents in each rack.

#### Parts of the 5-lane Bay

Both the *14-lane Bay* and the *5-lane Bay* dialogs display reagent information when a reagent rack or position is selected. In the *5-lane Bay* dialog, you can view the status of the racks, the reagents in each rack, and the remaining volumes of the reagents. You can also use this dialog to manually enter a reagent ID and volume for a position in a rack. Click the **5-lane Bay** button in the Machine Monitor to display this dialog. You can also access the *5-lane Bay* dialog from the *Resource Overview* window.

The tool tip (D, below) displays the *Use before* date. This indicates the first day which you can no longer use a particular reagent. The *Use before* date is the day immediately following the expiration date.

Rack Area		10	Reagent IDs	Reagent Prop	erties
1 2 3	4 5	Rack: S		Volume (ml)	Information
		1 E Reag1	492056110413643295	3.625	ABORH, Mono ctrl, 2011-02-26
		2 F Reag1	101140117221297392	6.841	7 ABORH, Anti-A, 2011-05-21
3 3		3 Reag1	203127112819822900	7.841	7 ABORH, Anti-B, 2011-05-08
		4 Reag1	504066117321629841	8.691	ABORH, Anti-D series 4, 2011-03-06
		5 Reag1	505127115710814010	8.725	? ABORH, Anti-D series 5, 2011-05
		6 Reag1	111289097180328537	8.365	? ABORH, A1-Cell, 2009-10-17
3 - 3	3 3	7 FReag1	113289097180504939	8.365	7 ABORH, B-Cell, 2009-10-17
Ē	B	8 Reag1	024282090450135192	6.610	mA Name: B-Cel 2001
	_	9 Reag1	230282090734074691	8.550	DA Use before: 2009-10-17 10
Done	Recal				Chibard the left: +

A: Lane button

B: Empty lane

C: Rack in lane

**D**: Tool tip – information about the reagent

The 5-lane Bay dialog is divided into three areas:

Rack

Reagent IDs

**Reagent Properties**
#### **Rack Area**

The Rack Area displays a graphical, aerial view of the 5-lane bay and the **Recall** button. The 5-lane bay graphic displays each of the 5 lanes with a corresponding button above. If a rack is present in a lane, the dialog displays a graphic of a rack in the lane. Click a rack to select it. A gray bar highlights the currently selected rack.

Use the lane selection buttons (1–5) to focus the barcode reader on the desired lane. When activated, the barcode scanner is focused on this lane and the instrument does not perform any operations on this lane. The colored bar in the button displays the status of the lane in accordance with the consistent color codes. The color displayed is the same as that of the lane LED on the instrument.

Click the **Recall** button to manually enter unreadable barcodes. If a reagent rack is loaded and one or more barcodes on that rack is/are unreadable, you can input the unread barcode. Refer to **Chapter 6** – **Instrument Testing Operation** for more information regarding how to rescan unread barcodes.

#### Reagent IDs Area

The Reagent IDs area displays details of the selected rack of the reagent loading bay graphic. The system displays the rack type in the Rack type field at the top of the Reagent IDs area, and reagent information in the list below.

The column contains tags that are used to select the type of material in the container.

The following options are available:



The software displays the reagent barcodes in the field of the Reagent IDs area. If the instrument is unable to read the reagent vial barcode, you can enter the ID manually. Refer to **Chapter 6 – Instrument Testing Operation** for information regarding the manual entry of unread reagent barcode data.

#### **Reagent Properties Area**

The Reagent Properties area lists the description of the reagent and also its remaining volume. By selecting the reagent description button, the tool-tip information about that reagent is displayed. The tool-tip lists the name, lot number, use before date, first on board, date and on-board time left, and QC status.

## Status Bar

### Using the Status Bar

### **Status Bar Buttons**

The Machine Monitor contains a Status Bar at the bottom of the screen. The system **Emergency Stop** button, in the middle of the status bar, stops all actions on the instrument at any time.



After you press the **Emergency Stop** button, the software requests confirmation. You must either press the **Yes** button to stop the system or the **No** button to continue NEO Iris operation. The NEO Iris will need to be initialized after the **Yes** button is pressed, and everything in process will be lost.



**Limitation**: Stopping the system with the **Emergency Stop** button interrupts instrument sampling.

The Status Bar also contains the three system status indicators that display the overall status for the instrument's connection to the PC, the level of the liquid waste in the common waste container, and the level of the available system liquid in the system liquid container.

Indicator	Description
¢	Connection of the PC to the instrument
Ъ	Level of the liquid waste collection
	Level of available system liquid

A status indicator has a gray background when the item has an acceptable normal state.



A status indicator has a yellow background when the item is reaching the point where action is recommended to either fill the system liquid because it is getting low or when the common waste container is nearing full and emptying the waste is recommended. The waste status indicator graphic also changes by showing the level of green liquid to have increased.



The system liquid indicator graphic also changes by showing the level of liquid in the bottles to have decreased to a very low level and a question mark is added to the graphic.



A status indicator has a red background when the item has reached a critical status point. For example: The NEO Iris is switched off and the PC that is still switched on cannot connect to the NEO Iris.



The common waste container is full and must be emptied.

The system liquid is insufficient and must be replenished.

### Words Indicating State of NEO Iris

The wording located between the Emergency Stop button and the three system status indicators describes the overall state of the NEO Iris. Wording options include *Not initialized, Initializing, Ready, Active, Halted, Starting Up* and *Login.* These states are described in the table below.

Wording	Description				
Not Initialized	Not Initialized: The NEO Iris is not initialized. No sample processing is permitted.				
Initializing	Initializing: The NEO Iris is in the process of initializing.				
0P Ready	<b>Ready</b> : The NEO Iris is initialized and is ready for use, but it is not currently active.				
	Note: Scheduled maintenance may still be due.				
Active	Active: The NEO Iris is initialized and is actively performing tasks such as processing samples, performing maintenance or scanning plate barcodes.				
Halted	Halted: The NEO Iris has been halted due to an unexpected event such as error recovery.				

**Starting Up**: This wording is only visible for a few seconds when the PC is initially powered off and is then powered on, but the database is not yet completely loaded. Once the database is completely loaded, a *Login* dialog is displayed and the **Login** wording is displayed, which replaces the **Starting Up** wording.

		T.

oain

Starting Up

The **Login** dialog and wording is no longer displayed after the operator logs in and the software display is completely established. Other wording is then displayed, depending on the status of the NEO Iris.



Note: This Login wording is only displayed under the scenario described here. When operators routinely log into and out of the NEO Iris software, while the NEO Iris software is displayed, then some wording other than Login is displayed, depending on the status of the NEO Iris.



**Note**: Another wording option is **Completing Run**. This wording option is only displayed when a software error occurs at the same time as when samples are being processed on the NEO Iris. That current run of sample processing is allowed to finish, however the starting of new processing runs is not permitted. The operator must shut down and then restart the software after that current run is finished.

### Log List

The *Log List* dialog displays a list of messages that have been issued since the last database clean up. You can also use this dialog to delete alarms or print a report of the message list.

### Viewing Messages

To view messages:

Step	Action
1	Press the <b>Log List</b> button (to the left of the <b>Emergency Stop</b> button) in the status bar.
	PC: LABTECH2 logged in with user level 1020000000     PC: Volanda logged out       PC: User administration: User LABTECH2' added by operator 'Yolanda'     PC: Beady

The	Log Lis	<i>t</i> dialog ap	pears.		
Sort	by: Date/Time	Searc	Log l	ist	
Туре	e Device	Number	Date/Time 🔻	Description	
	PC	00000708	2010-03-11, 17:34:48	VENDORSERVICE logged in with user level 1020000000	
O	PC	01000708	2010-03-11, 15:10:49	VendorService logged out	
0	PC	01000308	2010-03-11, 15:10:29	Machine is initialized	
	PC	00000308	2010-03-11, 15:08:21	Machine Initialization started	
	PC	0100070B	2010-03-11, 15:04:35	Error (Transport) occured - User VendorService pressed 'Abort run' for	Ξ
0	Transport	14030106	2010-03-11, 15:04:01	Y- Position Not Reached at position 08h params 01h while getting plate	
	PC	00000708	2010-03-11, 15:01:17	VENDORSERVICE logged in with user level 1020000000	
	PC	01000708	2010-03-10, 17:11:27	VendorService logged out	Ŧ
	PC	00000308	2010-03-10, 16:52:46	Machine Initialization started	
	PC	00000709	2010-03-10, 16:52:46	Event on connection 'Galileo Result Upload' - activated	T
	ок			Delete Detail F	rint

### Sections of the Log List

The *Log List* dialog is divided into three areas: Information area at the top Log list in the middle The buttons at the bottom

### **Information Area**

The information area at the top of the dialog comprises the **Sort by** and **Search key** fields. The Sort by field displays the name of the column that the event list is currently sorted by. To change the content of the Sort by entry, click on any of the field headers to sort the list by that field. The header of the field appears in the Sort by field. In the Search key field, you can enter search criteria for searching the log list.

### Log List

The Log List displays events and alarms that have occurred since the system was powered on. The system database stores entries until they are deleted. The log list contains five columns:

Type - displays an icon for the types of messages issued.

Device - displays the name of the system module that issued the message, for example, the PC.

Number - displays the message number, which is used to identify the entry.

Date/Time - displays the date and time the message was issued.

Description - displays a short description of the event that led to the message being issued.

Press the button at the head of each column to sort the list by that column. The Sort by field in the information area displays the name of the column that the list is currently sorted by.

### Message Types

The following message types are available:

Туре	Description
$\bigcirc$	General Message (event or warning)
0	Fatal Device Message
	Critical Device Message
8	Database
	Operating System
	Communication
8	Other

### Log List Navigation

On the right side of the dialog are seven navigation buttons that you can use to browse through the list. The following navigation buttons are available:

Button	Description
x	Go to the beginning of the log list.
\$	Go to the previous page of the log list.
	Go to the previous entry in the log list.
Ш	Select this button to mark multiple entries. Using this button, you can print or delete all marked entries as a group.
-	Go to the next entry in log list.
₹	Go to the next page of the log list.
¥	Go to the end of the log list.

### Log List Buttons

In addition to the common **OK** and **Print** buttons, the **Delete** and **Detail** buttons are also available. Press **Delete** to delete the entry selected in the log list. Press **Detail** to display the *Details* dialog.

		Detailed error description			ОК
Type: Number: Device:	Fatal Error 2C03020F Washer1	Date/Time:	2009-09-01, 14:20:59 0	_	
Descriptior Dispense 1=failed ro Top Prol Bottom Pro	: check failure : w HGFEDCBA bes 00010000 bes 00010000				
Recovery: Please ch	eck for proper priming, wa	ash head and dispense pump			
Initialise, tr	y again. If problem persist	is call support.			

The *Details* dialog displays more detailed information about the selected entry.

Press the **Print** button to display the *Export Selection* dialog.

### **Export Selection Dialog**

When one or more log list entries are selected and then the **Print** button is pressed, the **Export Selection** dialog is displayed. The *Selection mode* area allows filtering of the log list entries prior to export. One of the three (3) report options is *Excel Export*.

<ul> <li>Selected</li> </ul>	1 Items
C Range	By: Date/Time
◯ Age > 30 days	From:
C All	Τα
	Report
	Event Log, short
	Event Log, long



**Note**: This export process results in the selected log items being exported to the hard drive of the PC, and not to an external hard drive, memory stick or laboratory information system. The purpose of this localized export process is as a troubleshooting tool in the event that valuable information is saved in the event log lists. Such information can be then accessed remotely by Technical Support and used as part of the troubleshooting diagnostic effort.

Press the **OK** button to complete the selected export process or press the **Cancel** button to cancel the selected export process.

### **Details Dialog**

The *Details* dialog provides detailed information about the entry selected in the Log List dialog. Press the **Detail** button on the *Log List* dialog to display the *Details* dialog.

		Detailed error description			ОК
Type: Number:	Fatal Error 2C03020F	Date/Time:	2009-09-01, 14:20:59		
Device:	Washer1	Arm:	0 Sequence-ID:	0000	
1=failed ro Top Prol Bottom Pro	w HGFEDCBA bes 00010000 bes 00010000				
Recovery: Please chi	eck for proper priming, w	vash head and dispense pump			-
Initialise, tr	y again. If problem persi	ists call support.			

The detail screen displays the **Type** of message, the unique identifying **Number** of each error message, the **Device** involved, the **Date/Time** when the event occurred, and the **Arm**, **Block-ID**, and **Sequence-ID** indicating whether the pipettor was involved. The **Description** field contains a description of the event.

Press the **OK** button to exit the *Details* dialog.

# **Chapter 4: Security**

## In This Chapter

This chapter describes the procedures required to maintain NEO Iris users.

CHAPTER 4: SECURITY	4-1
Assigning Passwords and User Access Rights	4-2
Adding a User	4-3
Editing a User	4-10
Deleting a User	4-12
Default Access Rights	4-14
Changing a Password	4-16
Archive Configuration	4-18

## Assigning Passwords and User Access Rights

### Purpose

The NEO Iris system uses user names and passwords to track which user is currently operating the system and to restrict access to areas of the system by unauthorized operators.

Each operator has a unique user name and password with corresponding access rights. This section describes defining users with associated access rights, including adding a new operator and editing and deleting current operators.



**Note**: Even though the *Settings* dialog responds to touch-screen, Immucor recommends that you use the keyboard and trackball-mouse when using this dialog.

## Where to Do This

You administer user names, passwords, and access rights from the Security tab of the Settings dialog.



**Note**: Permission to use the *Security* tab is given by Immucor personnel to site operators on a case-by-case basis by configuring the site operator's user access rights.

## Adding a User

To add a new user:

Step	Action
1	Press the <b>Settings</b> button on the Main Menu Bar.
2	The system displays the <i>Settings</i> dialog.
	Press the <i>Security</i> tab.
3	The General Information tab is displayed as the default.
	Setup       Setup         Assay Definition       Profile Definition         Setup       Setup         Setup       OK         Setup       Password         Instrument Level : (000000       Op         Setup       Op         Cons       Op         Setup       Op
4	In the <i>Operator Name</i> field, type the new operator name.           Note         This field is not         case sensitive.

Step	Action			
5	Press the <b>Password</b> button to open the <i>Change password</i> dialog.			
	Change password         New Password         Enter Password :         Re-type Password :         OK			
6	Type the user password into the <i>Enter Password</i> field.			
	Note: This field is case sensitive.			
7	Re-type the user password into the <i>Re-type Password</i> field to confirm.			
	Note: This field is case sensitive.			
8	Press the <b>OK</b> button to accept the password entries. Both entries must be identical to proceed. Alternatively, press the <b>Cancel</b> button to cancel.			
	NoteNotePasswords must be between a minimum of five and a maximum of thirty alphanumeric characters in length.			
	Because manual entry of passwords in login dialogs is case sensitive, passwords must be entered exactly as set up for successful entry into the NEO Iris software.			
9	In the <i>Instrument Level</i> field, type the operator's instrument access level number. The access level number determines the access rights that are assigned to each operator. These access rights determine the functions of the software that the operator is authorized to access.			

Step	Action		
	<ul> <li>Note: The description and numeric assignments of the different access levels are as follows:</li> <li>Low (20000)</li> <li>Medium (30000)</li> <li>High (100000)</li> </ul>		
9 ( <u>cont</u> .)	<ul> <li>For access rights, see the section Default Access Rights in this chapter.</li> <li> i imitation: It is not possible to: <ul> <li>Assign to a new user higher user rights than the current user.</li> <li>Edit the user rights or change the password of a user with higher user r than the current user.</li> <li>Add a user who is already present on the user list.</li> </ul> It is important to consider the needs of each operator when assigning access rights. Setting access rights too low may prevent an operator from performing icb activities.</li></ul>		
	rights too high may permit an operator to gain access to portions of the instrument software he or she is unfamiliar with. Incorrectly assigning access rights can result in operator software changes that could negatively impact the performance of the NEO Iris.		
10	Press the <i>Menu Restrictions</i> tab.		



Step	Action			
13	The Assay Restrictions tab appears.			
	General Information       Menu Restrictions       Assay Restrictions       Language Settings         Enable individual Assay Restrictions       Enable individual Assay Restrictions       Select All         +       View Worklist Status       Select All         +       Add Worklist Request       Select None         +       View Sample Details       Select None			
	Real-time Sample Result Editing         +       2_Cell         +       3_Cell         +       4_Cell         +       4_Cell         +       ABORH         +       Ag_C RH2         +       Ag_c RH4         Vuser can perform function for assay         Please Note - These settings are specific to the 'Galileo_SQL' Database.			
	Two windows related to assay restrictions and result editing are accessible for configuring a user's assay ability. Each window has an associated check box. If the box is checked, then the contents of the associated window are editable with the correct security level. This activates the associated <b>Select All</b> and <b>Select None</b> buttons. If the box is not checked, then the contents of the associated window are not editable and the <b>Select All</b> and <b>Select None</b> buttons are inactivated. <i>Real-time Sample Result Editing</i> is used to configure the rights for real-time result editing of equivocal well results for certain assays. It is not possible for an operator to edit the final interpreted result of a sample. An example of not being able to edit an assay equivocal is shown below.			
	Real-time Sample Result Editing         CMV         CMV         Equivocal Wells         Press Select All to select all options. You can then selectively deny access to some specific         functions. Press the Select None button to remove all selections.			
14	Press the <i>Database Restrictions</i> tab.			

Step	Action		
15	The Database Restrictions tab appears.         Image: Definition Policy Po		
16	Press the <i>Language Settings</i> tab.		

Step	Action			
17	The Language Settings tab appears.			
18	Press the <b>OK</b> button to exit. All data entered is automatically saved.			

## Editing a User

To edit a user:

Step	Action		
1	Press the Settings button on the Main Menu Bar.		
	***		
2	The system displays the <i>Settings</i> dialog.		
	Press the <i>Security</i> tab.		
3	The General Information tab is displayed as the default.		
	Setup       Setup         Setup       OK         Setup       Devide Definition         Setup       Devide Definition     <		
	that user in the fields of the General Information box.		
4	<ul> <li>General Information</li> </ul>		
	Menu Restrictions		
	Assay Restrictions		
	Database Restrictions		
	Language Settings		

5	Press the <b>OK</b> button to exit. All data entered is automatically saved.

## Deleting a User

To delete a user:

Step	Action		
1	Press the <b>Settings</b> button on the Main Menu Bar.		
2	The system displays the Settings dialog.		
	Press the <i>Security</i> tab.		
3	The <i>General Information</i> tab is displayed as the default.		
	Setup       Setup         Image: Setup Device Device Setup Prote: Unit Meterinance Setup: US       Setup         Image: Setup Device Device Setup Prote: Unit Meterinance Setup: US       Setup         Image: Setup Device Device Setup Prote: Unit Meterinance Setup: US       Setup: UNIT Setup: Setup Prote: Unit Meterinance Setup: US         Image: Setup Device Device Setup Prote: Unit Meterinance Setup: US       Setup: Unit Meterinance Setup: Use Setup: Use Device Setup: US         Image: Setup Device Setup: Use Setup: Use Setup: Use Device Device Device Device Setup: Use Device		
4	Press the <b>Delete</b> button. The system displays the <i>User Warning</i> dialog that specifies the name of the user and asks "Are you sure?"		
5	Press the <b>No</b> button to cancel or press the <b>Yes</b> button to confirm. If you press the <b>Yes</b> button to confirm, the deleted user name disappears from the user list. A printer prompt appears with options to <b>Print</b> or <b>Cancel</b> . Printing will provide a log of the user being deleted.		
6	Press the <b>OK</b> button to exit.		

## Default Access Rights

The table below lists the Default access rights for the 20000 (Low), 30000 (Medium) and 100000 (High) User Levels.

The rights are as follows:

- ✓ Access Granted
- ⊗ Access Denied

Function	Low (20000)	Medium (30000)	High (100000)
Toolbar			
Login/Logout	~	~	✓
Initialize Instrument	√	✓	✓
Start Run	1	~	✓
Worklist	1	✓	✓
Plate Lists	1	~	✓
Results	√	✓	✓
Maintenance	√	~	✓
Utilities	8	✓	✓
Settings	√	~	✓
Help	1	~	√
Exit	√	~	√
Module Dialogs			
Washer	~	~	✓
Incubator	√	~	✓
Reader	1	~	✓
Loading Tower	~	~	✓
14-lane and 5-lane Bays	~	~	✓
14-lane Bay barcode configuration	8	✓	✓

Function	Low (20000)	Medium (30000)	High (100000)
Start Run			
Sample Loading	✓	✓	✓
Resource Loading	✓	✓	✓
Plate List			
Processing Steps	✓	✓	√
Delete	8	8	8
Detail View	✓	√	√
Print	8	8	√
Settings			
Assay Definition	8	8	8
Profile Definition	8	8	✓
Setup Barcodes	8	√	✓
Setup Device	8	8	8
Setup Microplates	8	8	8
Users	8	8	8
Maintenance	8	8	8
Service	8	8	8
Plate Loading Tower			
Strip Selection	1	✓	✓
Scan Plate	✓	✓	✓
Assay Selection	✓	✓	✓
Plate Expiry Entry	✓	✓	✓
Plate Expiry Removal	8	8	✓
Processing Steps	✓	✓	✓

## Changing a Password

You can change your password without having full access to the security configuration as long as your access has been set up by another operator with appropriate access rights of their own. Your access rights must be configured to *Change Password*.



To change your password, without rights to changing access rights:

Step	Action			
1	Press the <b>Settings</b> button on the Main Menu Bar.			
	The Security tab will open as the default and the Change password dialog automatically			
	displays. No additional operator set up options will be available.			
	Current Password			
	New Password Enter Password : Re-type Password :			
	OK Cancel			
2	Enter your current password into the <b>Current Password</b> field followed by two (2) identical			
	entries in the two (2) fields in the <i>New Password</i> area (Enter Password and Re-type			
	Password).			

Press the Cancel button to cancel your entries, or press the OK button to register your new password.
 The Change password window will close and the blank security tab will be displayed indicating that no options are available to you.
 <u>Note</u>: You will have to use your new password when you log on next time.

## Archive Configuration

This section describes how to configure the archive process and view archive logs.



**Note**: The purpose of describing the archive configuration is to inform you of the configurations available. An Immucor employee with sufficiently high security access rights is required to configure the available options. You will not be able to access these configurations.

## **Configuring the Archive Process**

To configure the Archive process:

Step	Action				
1	The <b>Utilities</b> dialog allows you to view event logs and statistics, as well as archive and print reports. Select the <b>Utilities</b> button to access the <b>Archive</b> tab.				
	Event log       Archive       Statistics       Reports       Utilities         Drive space:       Final archive location       Progress       Current action :       Total :         Drive space:       Free:       702,82 MB       Total :       Total :       Total :         Drive space:       Free:       72,92 MB       Total :       Total :       Total :         Drive space:       Free:       51,97 GB       Total :       Total :       Total :         Drive space:       Free:       51,97 GB       Total :       Total :       Total :         Drive space:       Free:       51,97 GB       Total :       Total :       Total :         Estimates       Space:       Free:       54,23 GB       Total :       Total :       Total :         Status       Oldest existing data:       2009/08/14       Last archive date:       2009/08/13       Total :       Total :         Norting database:       GalleoSQL       Archive       Option of y       Onfigure archive arc	0% 0% 601			

Step	Action
2	The <b>Configure archive</b> button is only available on the <b>Archive</b> tab if the Immucor employee security configuration allows access.
	Configure
	Archive Perform Archive Configure Archive
	Security access permissions for archiving configuration
	Note:The Action and Archive boxes are not normally accessible on the Archive tab, but show what has been selected on the configuration screen. Access to these boxes can be activated from the General options tab as described later.
3	The Archive settings dialog is displayed when the Configure archive button is selected.
	Archive settings       ? ×         General options       Schedule next archive       Access specific       SQL Server specific         Archive options       Archive options       Image: Transmission of days kept in the working database : 7 *       Image: Transmission of days kept in the working database : 7 *         Archive actions       Image: Transmission of days kept in the working database : 7 *       Image: Transmission of days kept in the working database : 7 *         Image: Copy records to archive database       Image: Transmission of Data and Images : Transmission of Data only       Image: Transmission of Data only         Image: Allow action selection before archiving       Image: Transmission of Data only       Image: Transmission of Data only         Image: Copy records to archive target path : Image: Transmission of Data only       Image: Transmission of Data only       Image: Transmission of Data only         Image: Distribution of Data on the Data

The Archive settings dialog has four tabs:

- General options
- Schedule next archive
- Access specific
- SQL Server specific

### **General options tab**

A number of configurations can be set up through the general options tab. The number of day's worth of data that should remain in the working database after an archive, assuming records have been moved (copy and deleted), can be set using the **Number of days kept in the working database:** numeric selection box. Seven (**7**) days is the maximum number of day's worth of data that can be kept. Zero (**0**) is the minimum number of day's worth of data that can be kept. Zero (**0**) is the minimum number of day's worth of data that can be kept. You can select whether to copy or move (copy and delete) records to the archive database. You can configure whether the move/copy data to the archive database can be selected on the main archive dialog using the check box for **Allow action selection before archiving**. If selected, the **Action** and **Archive** boxes are displayed on the main archive window.

Action ———	Archive
Move	O Data and images
С Сору	🔿 Data only
C Delete	

Action box

The **Move** action is defined as copy and delete. The **Copy** action is defined as copy without delete. The **Delete** action is defined as delete without copy.

Archive box



**Note:** The **Delete** option of the **Action** box is only active if a copy action has been performed.

You can select whether to archive **Data and images** or **Data only**. If the selection is **Data only** and the action is to move (archive) data, any associated image files within the corresponding date range will be deleted.

You can define the final archive target path, e.g. another area of the hard disk, network or removable media. The correct target is E:W.

The Immucor employee can accept the required configurations by selecting the **OK** button or reject any changes by selecting the **Cancel** button.

### Schedule next archive tab

The **Schedule next archive** tab allows the archive interval to be defined. The default and maximum allowable schedule time is **7** days. Archiving can be performed more frequently than 7 days if needed or preferred. The allowable numeric range for the number of days selection is between one (1) and seven (7) days inclusive. Zero (0) days and greater than seven (7) days cannot be selected.

General options	Schedule next archive
O Days	7 🛖

Schedule next archive tab

The instrument will not prevent you from running assays if this scheduled time is exceeded. However, there will be a warning flag generated on the Archive tab within the Utilities dialog and also on the Utilities button.



Refer to **Chapter 10 – Maintaining the NEO Iris** for more information regarding the incorporated archive warning system.

The Immucor employee can accept the required configurations by selecting the **OK** button or reject any changes by selecting the **Cancel** button.

### Access specific tab

The **Access specific** tab is not needed for routine archive configuration. The entries in this tab should not be modified without consulting Immucor.

The Access specific tab allows definition of a local archive path. The default value is

**D:**#MPA#data#Aurora#Archive. The check box is used to configure **Compact and repair databases** when archiving completes.

General options   Schedule next archive	Access specific	SQL Server specific
Local archive target path :		
Compact and repair databases when archiving completes		

#### Access specific tab

The Immucor employee can accept the required configurations by selecting the **OK** button or reject any changes by selecting the **Cancel** button.

### SQL Server specific tab

The **SQL Server specific** tab is not needed for routine archive configuration. The entries in this tab should not be modified without consulting Immucor.

#### **Chapter 4: Security**

The **SQL Server specific** tab allows definition of a server name for archiving. The drop down list finds and lists any servers present. Once the server is selected, a **Verify server** button will become active with a warning symbol. When selected, the software will verify that the server name entered exists and that it is configured for remote access. If successful, the warning symbol will then be removed.

General options Schedule next archive Access specific SQL Server specific		
Archive database settings		
Server : svrsource	<b>•</b>	
	<u> </u> Verify server	

SQL Server specific tab

A check box for **Shrink/truncate databases when archiving completes** is located on this SQL Server specific tab.

The Immucor employee can accept the required configurations by selecting the **OK** button or reject any changes by selecting the **Cancel** button.

## **Viewing the Archive Logs**

A list of archive activity can be viewed using the archive logs. The **Archive Logs** are accessed by pressing the **Archive logs** button on the main *Archive* dialog.



Archive logs button

The Archive Logs dialog has three columns:

- Date
- Operator
- Action

The *Date* column documents the date when the archive activity took place. The *Operator* column documents the identification of the operator who performed the archive. The *Action* column describes the details of the archive activity.

ICHIVE LU	15		
Date	Operator	Action	_
2010/02/0	2 Yolanda	Archiving Version 2 Type: Data only Action: Copy data Keep: 7 days data (<=2010/01/26) Target Name: 20090929_20100126 Server: Galileo-PC	
2010/01/2	7 Yolanda	Final: E:\ Archiving Version 2 Type: Data only Action: Copy data Keep: 7 days data (<=2010/01/20) Target Name: 20090929 20100120	
			Close

When large quantities of entries are listed, horizontal and vertical scroll bars can be used to navigate the **Archive Logs** dialog. Press the **Close** button of the **Archive Logs** dialog to exit the dialog.

# **Chapter 5: Instrument Start-Up**

## In This Chapter

This chapter describes how to start up the instrument.

CHAPTER 5:	INSTRUMENT START-UP	5-1
Starting Up	)	5-2
Logging In	and Initialization	5-3
Auto Logo	ff	5-5

## Starting Up

### Purpose

To use the instrument, you must turn on the NEO Iris PC and the NEO Iris module.



**Note**: The centrifuge module should already be switched on unless it was switched off during some troubleshooting process.

These steps should be followed when the NEO Iris instrument is needed for test processing or to perform the daily shutdown/startup of the instrument.

## Procedure

To start the NEO Iris:

Step	Action
1	Turn ON the computer.
	Note: Open the right door of the cabinet to locate the computer.
2	Turn ON the monitor.
3	Turn ON the NEO Iris instrument.
	<b><u>Note</u></b> : The power switch is located next to the monitor pole.
	Note: Although the centrifuge module is not typically powered off, it may be necessary to turn ON if the module was turned off during the shutdown procedure. The power switch is located at the rear of the module.
4	The Camera Reader lamp and instrument overhead lights will be illuminated when the NEO Iris is powered ON.
	Limitation: The instrument must be switched on at least thirty (30) minutes prior to the first plate read to allow the reader lamp to warm up. Reading of plates prior to completion of this warming period can cause incorrect negative reading of weakly positive reactions.
5	Log into the NEO Iris software.
	Refer to Logging In and Initialization published in this chapter for detailed information.

## Logging In and Initialization

### Purpose

To access the instrument software, you must enter a valid operator ID and password. For more information about operator IDs and passwords, refer to **Assigning Passwords and User Access Rights** in Chapter 4 - Security.

**Note**: The instrument software is password protected to prevent unauthorized operator access to the computer operating system and also prevents configuration changes by unauthorized personnel.

### Procedure

Step	Action
1	The Login dialog will be displayed either automatically at startup or by selecting the Login button on the Main Menu Bar.          NEO Iris - Login         Please enter your name and password         Please enter your name and password         Password :         Start Action         Restore Databases from Backup         Use DMS Archive         OK       Shutdown
2	In the Name field, type your operator name.           Mote:         This field is not case sensitive.
3	In the Password field, type your password.           Mote         This field is case sensitive.

Step	Action
4	If you have sufficient user rights, you can select the <b>Use DMS Archive</b> check box.
	Limitation: The Use DMS Archive check box control option allows you to view archived results that are saved on a disk or another Archive medium using the DMS. You can only view archived results through the DMS with the Archive medium attached to the computer or by inserting the disk into the drive of the computer.
5	Click <b>OK</b> to complete the <b>Login</b> process.
	The system automatically performs an initialization routine that includes removing any plates that may have been left on the instrument and returns them to the Plate Loading Tower.
	<b>Attention</b> : If the instrument displays a Centrifuge: Rotor error message during initialization, the centrifuge module is not powered on.
	Limitation: The COP: Serial Buffer Deleted error message may appear during initialization. In this case click the <b>OK</b> button. However, if this error occurs in other situations it is recommended to complete the current run(s) but do not start any new plates. Initialize the instrument using the Initialize button before starting new plates.
6	The instrument is now ready to perform system maintenance and sample processing.
### Auto Logoff

#### Purpose

If the automatic Logoff function is configured, the system will activate a locking screen after a certain period of inactivity. You must enter a valid operator ID and password to access again the instrument software.



**Note**: The optional Auto Logoff function is only available with Instrument Control Software version 1.8 and higher.

#### Procedure

Step	Action
1	The AutoLogoff dialog will be displayed automatically after the configured duration of the inactivity period will elapse.          AutoLogoff         Please enter your name and password         Name :         Password :
2	In the Name field, type your operator name.           Mote         This field is not         case sensitive.
3	In the Password field, type your password.           Mote         This field is case sensitive.

Step	Action
4	Limitation: Once the Auto Logoff has occurred, only users with equal or higher access rights than the user that was logged in when the Auto Logoff was triggered are able to access the system.
5	Click <b>OK</b> to unlock the system.
	The instrument software is again accessible.

# Chapter 6: Instrument Testing Operation

#### In This Chapter

This chapter describes the detailed operating procedures required to process samples.

CHAPTER 6: INSTRUMENT TESTING OPERATION	6-1
Using the Start Run Assistant	6-2
Loading Samples	
Downloading Requests from LIS	6-14
Completing the Sample Loading Process	6-16
Loading Reagents and Controls	6-21
Loading Plates	6-27
Starting Processing	6-34
Continuous Loading During Operation	6-36

### Using the Start Run Assistant

#### **Before You Begin**

You must prepare all of the necessary reagents and samples according to the reagent package insert requirements before running assays. The *Start Run Assistant* dialog provides an intuitive guide to starting assay runs.

#### About the Start Run Assistant

Click the Start Run Assistant button on the Main Menu Bar.



The instrument displays the Start Run Assistant dialog.

Start Run Assistant	
1. If specimen tubes have not yet been loaded please go to "Load Samples" .	Load Samples
2. Query the host system for test orders.	Download Requests
3. Complete loading (micro plates, reagents, controls, disposable tips, and wash buffers).	Load Resources
Cancel	

#### Start Run Assistant Steps

The following steps must be completed:

- 1. Load Samples
- 2. Download Requests from LIS (optional)
- 3. Complete the Sample Loading Process
- 4. Load Reagents and Controls
- 5. Load Plates
- 6. Start Processing

### Loading Samples

#### Sample Requirements



Refer to **Attachment 1 for Neo Iris Operator Manual** for details regarding sample requirements.

#### Procedure

To load samples:

Step	Action
1	Click the <b>Load Samples</b> button on the <i>Start Run Assistant</i> dialog. The instrument displays the <i>14-lane Bay</i> dialog.
	Prefix + Barcode Setting:       Test Selection         2 3 4 5 6 7 8 9 10 11 12 13 14       Prefix + Barcode Setting:       2 Cell         4 Dell       Ab_D         4 Dell       Ab_D         4 Dell       Ab_D         ABORH       ABORH         Brev       S TAT         Next       Profile
2	Ensure that the sample tubes are seated in the appropriate rack and the barcode labels can be seen through the gap on the right of the sample rack.
	Note: Correct position requires that the tubes fit securely, with the barcode positioned between 20 mm and 105 mm (0.8 and 4.1 inches) from the bottom of the tube.

Step	Action
3	Slide the sample racks into the sample loading bay one at a time. Use the lane with the continuous green indicator LED.

Step	Actio	n
4	Verify that all the barcodes have been scanned su appear in the fields of the Identifiers area of the Rack: B Prefix + Barcode Settings 1 Sample 12202 2 Sample 12205 3 Sample 12204 4 Sample 12203 5 Sample 12206 6 Sample 12206 6 Sample 2345677 9 Sample 2345677 9 Sample 2345678 10 Sample 456789 11 Sample 456789 1	uccessfully. Successfully scanned barcodes <i>14-lane Bay</i> dialog.
	If	Then
	All barcodes have been scanned successfully	Load further racks of samples as necessary.
	A barcode has not been scanned successfully	Follow the procedure to <b>Rescan Unread</b> Sample Barcodes.

#### **Barcode Settings**

The *Barcode Settings* dialog can be used if sample barcodes require any data conversion or manipulations on a rack-by-rack basis. This dialog is accessed by pressing the **Prefix + Barcode Settings** button located top center of the sample loading dialog.



- Sample tubes requiring barcode conversion should be loaded into racks with tubes of like barcode properties
- Configurations defined via this dialog will remain active until the 14-lane bay is exited or the settings are modified again

- If racks that require different barcode conversions, the Prefix + Barcode Settings dialog must be opened to modify the current conversion scheme before loading a rack with a different conversion requirement
- To cancel the current conversion, the dialog can either be opened and conversion de-selected or exit the 14-lane bay to revert back to any default barcode settings.

The Barcode Settings dialog has two tabs (Donor and Sample).

code Settings	Barcode Settings
onor Sample	Donor Sample
Mask:	Mask:
Prefix:	
Codabar Conversion:	Codabar Conversion:
C ABC	C ABC
C Custom	C Custom
ISBT 128 Conversion	ISBT 128 Conversion
Ok Cancel	Ok Cance

Barcode Settings>Donor tab

Barcode Settings>Sample tab

The Donor and Sample tabs have the same features for being able to set a barcode mask, converting Codabar (ABC or Custom) and ISBT 128 conversion. The Donor tab has the additional feature of being able to assign a designated prefix to all samples on a given donor rack. Separate prefixes cannot be assigned to different samples in the same donor rack. The maximum number of characters that can be entered into the Prefix field is five. The tabs will populate with the pre-configured settings that were previously inputted into the *Setup Barcodes* tab of the *Setup* screen, however this default can be changed by you on a rack-by-rack basis.

If a donor prefix is configured, it will be displayed as replacement text on the *14-lane Bay* dialog *Prefix + Barcode Settings* button. For example, if a donor prefix of *12345* is used, then it will display, as shown below, as replacement text on the button.



Press the **Cancel** button of the *Barcode Settings* dialog to cancel any selected options and close the dialog. Press the **OK** button of the *Barcode Settings* dialog to accept any selected options and close the dialog. Refer to **Chapter 3 – System Software Navigation** for details regarding customized Codabar conversion and masking.

### **Rescanning Unread Sample Barcodes**

To rescan unread sample barcodes:

Step	Action	
1	Remove the rack with the unread barcode and ch	eck that the barcode is correctly aligned.
2	Slide the rack back into the loading bay. Use the LED.	lane with the continuous green indicator
3	Verify that all the barcodes have been scanned su appear in the fields of the Identifiers area of the Z Rack Prefix + Barcode Settings 1 Sample 1200 2 Sample 12205 3 Sample 12205 6 Sample 12206 6 Sample Colin123 7 Sample Beatrix234 8 Sample 2345678 10 Sample 3456789 11 Sample 4567891 12 5 13 5	ccessfully. Successfully scanned barcodes
	If	Then
	All barcodes have been scanned successfully	Assign required assays.
	A barcode has not been scanned successfully	Go to step 4.

4	a. Remove the rack and select the <b>Recall</b> button to display the sample IDs on the screen for the last rack removed.
	b. Hand-scan or type any missing / unread barcodes into the relevant field (double data entry is required)
	c. Reinsert the rack in the lane with the continuous green LED
	Note: You can make assay assignments at this time if you are assigning assays manually (within the 14-lane Bay dialog). If you have requested the assays through a LIS, close the 14-lane Bay dialog and proceed to <b>Downloading</b> Requests from LIS.
	If the assays have been requested using the Worklist facility, then the assays will be automatically assigned when the samples are loaded.

### Manually Selecting the Assays Required

To manually select the assays required:

Step	Action
1	Select the rack to receive assay assignments within the 14-lane Bay dialog.
2	Select a sample ID to choose a specific sample. The field is highlighted in blue. For STAT samples, proceed to step 3. For non-STAT samples, proceed to step 4.



Step	Action		
4	For non-STAT samples, press an assay button or assay profile button in the Test Selection area to select an assay or assay profile.		
	Note: If a crossmatch assay is selected, then an additional input window is displayed (titled with the name of the assay) after the assignment of that assay to a specific sample. The ID of the specifically selected sample is displayed just below the name of the assay.		
	IoB_XM         Beatix234         Selected Donors:         1234567         2345678         3456789         4567891         Delete         New Donor:         No X-match         OK       Cancel		
	This window allows the operator to input the specific donor unit segment IDs that are to be crossmatched with that specific sample.		
	The operator must hand scan or manually enter the donor ID(s) information into the <b>New</b> <b>Donor</b> field:		
	a. Hand-scanning will automatically add the Donor ID to the Selected Donors dialog; manual entry of id must be followed by pressing the <b>Enter</b> key on the keyboard to add the Donor ID to the Selected Donors dialog.		
	b. After all donors are in the Selected Donors box, select <b>OK</b> and then load the donor rack onto the instrument and move onto the next step.		
	Note:Nete: The Prefix field is populated by the configured value from the SetupBarcodes tab of the Setup screen. As barcode data is entered into the SelectedDonors field, the prefix is automatically assigned to the ID in the list.		

Step	Action
	Use the <b>Delete</b> button of the input window to remove specifically highlighted donors ID(s) from the <i>Selected Donors</i> list. Press the <b>Cancel</b> button of the input window to cancel the current donor assignments. Press the <b>OK</b> button of the input window to accept the donor assignments. If crossmatches have been previously assigned and accepted by pressing <b>OK</b> , and you wish to subsequently remove them, you can reenter the crossmatch input window and press the <b>No X-match</b> button to inactivate the assignments.
	Note: You must press the <b>OK</b> button to close the dialog before loading the donor rack into the bay, otherwise the tests will not be ordered correctly.
4 ( <i>cont.</i> )	Other assays, such as red blood cell antigen screening assays, utilize equivalent data entry windows with the same kind of functionality. The data entry window is titled with the name of the assay. Examples are shown below:          mAgScrII       mAgScrII       PAgScrAHG         Beatix345       Selected Donors:       Prefix         1234567       Prefix       Prefix
5	Repeat steps 1–4 for all samples requiring that assays be manually assigned.         Image: To apply the same assay selection to all samples in a rack, place the curser in the first position and then select the All button. The first position can be empty.         If you select a position other than the first position and click All, the assay will NOT be ordered for all samples.
	Refer to <b>Attachment 1 for Neo Iris Operator Manual</b> for a full list of assays and their respective assay button abbreviations. Profile buttons are site-customized buttons that can be built to combine two or more assays. Profiles are a group of assays that can be selected together and therefore make assay selection easier.

Action
You must click <b>Done</b> on the <i>14-lane Bay</i> dialog when sample loading and assay assignment is finished. The system redisplays the <i>Start Run Assistant</i> dialog.
Start Run Assistant
1. If specimen tubes have not yet been loaded please go to "Load Samples".
2. Query the host system for test orders.
3. Complete loading (micro plates, reagents, controls, disposable tips, and wash buffers).
Cancel
Note: The selected check box to the left of the Load Samples button indicates

# Downloading Requests from LIS

#### Purpose

You can download all test requests for samples from the LIS if a suitable interface is active.

The configuration of the LIS interface may require manual query for assay assignments using the procedure below or the assay assignments may be received by the instrument in real-time and therefore immediately appear when the sample is loaded on the instrument.



**Note**: You must validate the interface communication between the LIS and the instrument software prior to use. It is your responsibility to perform this validation. For more information on the configuration requirements, refer to **Attachment 2 for Neo Iris Operator Manual**.

#### Procedure

To download Assay Selections from the LIS:

Step	Action
1	Click <b>Download Requests</b> on the <i>Start Run Assistant</i> dialog.
2	The instrument queries the LIS for assay requests for the samples loaded and then downloads these requests from the LIS.
3	When it has finished, the system redisplays the <i>Start Run Assistant</i> dialog.

4 If the instrument is currently running assays, Open and close the *Loading Bay* dialog (via **Load Samples** of the *Start Run Assistant*) **AFTER** the LIS worklist orders have been received but **BEFORE** the *Resource Overview* window is opened. If no assays are currently running, skip to *Completing the Sample Loading Process* section.



<u>Note</u>: The selected check box to the left of the **Download Requests** button indicates this stage has been completed.



**Limitation**: When a sample is actively being tested by an assay and the same assay is requested again for that same sample from the LIS, a second request is added. This will cause the sample to appear in the *Resource Overview* window for a duplicate test. To avoid this, open and close the *Loading Bay* dialog (via **Load Samples** of the *Start Run Assistant*) **AFTER** the LIS worklist orders have been received but **BEFORE** the *Resource Overview* window is opened.

# Completing the Sample Loading Process

#### Procedure

To complete the sample loading process:

Step	Action
1	Click <b>Load Resources</b> on the <i>Start Run Assistant</i> dialog. The system displays the <i>Resource Overview</i> window.
	Resource Dverview     Plates     Reagents     Controls     Donors     Washbuffer     Pipettor     Incubator       DAT     1/1
	Pool_Cell       1/1         ●       Reflex&B0       1/1         ●       Reflex&D       1/1         ●       Reflex       ■         ●       Reflex       ■         ●       Reflex       ■         ●       Reflex       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■         ●       ■       ■      <
	PC: Scanner initialized PC: Scanner initialized PC: Scanner initialized Sample Loading Bay: Scanner time-out
	The <i>Resource Overview</i> window presents a list of plates (with the assay name and number of samples on each) that need to be processed to fulfill the assays ordered for the samples currently loaded onto the instrument.
	For information regarding STAT test processing, see below. To continue to the next step in this procedure, see the next page.

#### STAT Test Processing

The Resource Overview window preferentially displays STAT test orders before routine test orders. Existing routine test orders do not display until STAT tests start. After STAT tests start, routine test orders reappear in the window and their tests can be started.

- When other normal priority samples on the instrument require the same test as the STAT sample, the instrument uses the remaining positions on the plate for those samples.
- The Resource Overview window's *Samples/Strips* column displays:

o X(+Y)/Z where:

- X is the number of STAT tests ordered
- Y is the number of normal priority tests ordered
- Z is the number of strips required or loaded to run the STAT tests
- After you select a Resources Overview window line, the software checks available
  resources (including plates) and if STAT samples are included, it recalculates how it will
  use them. The *Samples/Strips* display changes, based on strip availability and resources
  needed. The number of normal priority samples displayed (Y) can decrease, based on the
  remaining capacity of the strip wells available after STAT samples are accommodated.

Samples/Strips Display Examples After STAT Tests Ordered	Number of Strips Loaded	Status Examples
ResourcesAssay NameSamples/StripsPlat3_Cell $1(+4)/1$ STAT $X(+Y)/Z = 1(+4)/1$ X = 1 STATY = 4 normal priority samplesZ = 1 strip required for the 1 STAT	None	Assay cannot be started because no strips are available. If one strip is loaded, it will be used to run the STAT sample and as many of the normal priority samples as can be accommodated, based on the remaining strip capacity.
ResourcesAssay NameSamples/StripsPlaAg_C RH2 $1/1$ STAT $X(+Y)/Z = 1(+0)/1 = 1/1$ X = 1 STATY = no normal priority samplesZ = 1 strip required for the 1 STAT	None	Assay cannot be started because no strips are available. When one strip is loaded, it will be used to run the STAT sample only.
Resources Assay Name Samples/Strips Plate 3_Cell 1(+1)/1 STAT X(+Y)/Z = 1(+1)/1 X = 1 STAT	1 strip	Assay can be started. The one strip will be used to run the STAT sample and include one normal priority sample.

Y = 1 normal priority sample	
Z = 1 strip loaded for the 1 STAT	

### **Procedure (continued)**

Step	Action
2	To select a plate to be run, click anywhere in the line of the screen where the plate is displayed. The system calculates and displays the status of the consumables (plates, reagents, controls, and wash buffer) for that plate.
	If the system displays a green check mark ( $\ll$ ) in the column of a consumable, this indicates that sufficient resource is available to perform the assays on the selected plate.
	Image: Section ABO         1/1         Image: Section ABO         1/1         Image: Section ABO         I
	If the system displays a red exclamation symbol ( <b>0</b> ) in one of the columns for a consumable, you must load that consumable.
3	Open the <i>Loading</i> dialog for the missing resources directly from the <i>Resource Overview</i> window by clicking on the appropriate button above the column with the exclamation mark:
	Plates
	Reagents
	Controls
	Washbuffer
	Pipettor
	Incubator
4	Resolve the resource issue using the missing list provided by the software.

Step	Action
5	Close the resource dialog. The system redisplays the <i>Resource Overview</i> window.
	Resource Overview       Plates       Reagents       Controls       Donors       Washbuffer       Pipettor       Incubator         DAT       1/1
	PC: Scanner initialized PC: Scanner initialized
	Sample Loading Bay: Scenner time-out You cannot adjust the number of strips on a plate in the plate loading dialog, if you have already selected that plate through the <i>Resource Overview</i> window. You must deselect the plate in the <i>Resource Overview</i> window before you can change the number of strips on the plate.
	Note: You should not remove sample or reagent vial racks while the <i>Resource</i> <i>Overview</i> window is currently being displayed. If a rack is removed, a message will

appear and Cancel must be selected to start the loading process again.

# Loading Reagents and Controls

#### Procedure

To load Reagents and Controls:

Step	Action
1	Load the reagents and controls required for the planned assays into the instrument reagent racks. Reagent rack details are listed in <b>Appendix A – Preparing the NEO Iris for First Use</b> .
	Note: The vials of cellular reagents should be gently mixed by hand, before removing the dropper assembly or cap, to completely resuspend the cellular contents (prior to adding the stirball in step 3 below).

Step	Action
2	Check the position of the vials in the reagent rack to make sure that the NEO Iris designated vial barcodes on the vial label are facing towards the opening in the right hand side of the rack.
	<b><u>Limitation</u></b> : Excluding QC assays, some assays will not be processed under the following conditions:
	A corQC EXTEND Standard vial loaded in the 5-lane bay will be recognized; however this reagent will be listed as missing in the <i>Resource Overview</i> window due to the risk of a probe arm crash.
	When running a DAT assay, a DAT Positive Control Cell vial loaded in the 5-lane bay will be recognized; however this reagent will be listed as missing in the <i>Resource Overview</i> window due to the risk of a probe arm crash.
	Under both of these scenarios, the assays will not be able to be scheduled until the necessary reagent is loaded into the 14-lane bay.
	Warning: Inspect all reagents and controls for the presence of foam before placing on the instrument. Do not vigorously agitate blood grouping anti-sera or controls. Shaking will produce foam in the vial that can cause the Liquid Level Detection (LLD) feature of the pipetting system to aspirate foam and/or air rather than reagent. This will produce incorrect results or an error.
	Warning:Before placing reagents on the instrument, you must remove the bottle caps. You are advised to remove and discard the dropper by pulling the dropper from the bulb. When you remove the reagents from the instrument for storage, you must place the caps back on the bottles. To avoid cross contamination of reagents, it is important that you place the caps on the correct bottles. Mixing caps can result in erroneous test results. Immucor-approved disposable vial caps can alternatively be used instead of replacing the saved original vial caps. Contact
	Technical Support for information regarding these caps.

Step	Action
3	Add one stirball into each vial containing cellular reagents, such as Referencells, Capture Indicator Cells, and quality control cells.
	A stirball is not required if one has already previously been added to a given vial prior to previous usage.
	<b>Note</b> : The stirballs are used to keep cellular reagents in suspension during testing.
	Warning: If you do not add the stirballs to the cell suspensions, the results may be invalid or incorrect. Do not touch the stirballs. You should add them directly to the cellular reagent vials using the dispenser provided. Contamination and neutralization of cellular reagents can occur if the stirballs are touched.
	You must only add one stirball per vial of cellular reagent. Do not add more than one stirball per vial.
4	Depending on the reagent rack size being used, open the <i>14-lane Bay</i> by selecting the Controls button at the top of the <i>Resource</i> Overview window or open the <i>5-lane Bay</i> dialog selecting the <b>Reagents</b> button at the top of the <i>Resource Overview</i> window.
	B: Empty lane

Step	Action
5	Load the reagent racks into appropriate slots following the indicator lights. You can introduce racks into the loading bay when the indicator light is continuous green. You can place non-cellular reagents and controls in any lane.
	Note: You should place racks that are accessed using the single probe arm to the right hand side of the loading bay; for example, a reagent rack. You should place a rack that is accessed only by the 4-probe arm on the left-hand side of the bay; for example, a sample rack. This practice optimizes the instrument's operational efficiency.
	You must place racks that contain cellular reagents in the lanes that have stirrer mechanisms (two positions on the right-hand side of the 14-lane Bay).
6	You can change the lane that is activated for insertion of a rack by clicking one of the indicator buttons at the top of the lane as represented on the <i>14-lane</i> or <i>5-lane Bay</i> dialog.
	The button above the selected lane changes to continuous green color, indicating that the barcode scanner is activated for this lane.

Step	Action
7	Wait for the indicator light (LED) at the front of the lane itself to change to continuous green before inserting a rack. This may take a few seconds.
	Limitation: When you insert a rack into the activated lane, the system scans the barcodes. The system displays the barcode data and interpretation (reagent name, lot number, expiration date) on the right-hand side of the on-screen dialog. The system also displays the volume remaining in each container. In the case of Immucor lot numbered reagent vials, the volume is based on the full volume of the container (for new vials) or the remaining volume (if the vial was previously used on the instrument). If the actual volume is less than that displayed (for example, if the reagent vial has also been used for manual testing), you may enter a reduced volume. It is not possible to enter an increased volume for the bottle.
	Warning: If you are using two or more instruments, then the specific reagent vials for each instrument must be dedicated for use on that single instrument to ensure correct reagent volume tracking. If the actual reagent volume (less than the software numeric volume) is not sufficient for the number of tests scheduled, the instrument will produce invalid results and samples will need to be rescheduled for testing.
8	Verify that the instrument has read all barcodes. If a barcode is not read, you may remove the rack, press the <b>Recall</b> button, hand scan the necessary barcode into the correct field (double entry required), and reinsert the rack.
9	If the entries are inconsistent, the system displays a <b>Warning</b> message and you must re-enter the duplicate entries.
	Note: This method uses the <b>Recall</b> button. After removing the reagent rack, click the <b>Recall</b> button to display the reagent IDs on the screen of the last set of barcodes read. You can then scan any missing unread barcodes into the relevant field with the hand held scanner.
	Note: If an assay consumes large volumes of a reagent (for example, Capture LISS), you can load several vials of the same reagent. The instrument will automatically use another vial when the first is empty.

Step	Action
10	When you have finished loading reagents, click <b>Done</b> to exit the <i>Loading Bay</i> dialog.
	Limitation: If a reagent expires during a testing process, the results will not be flagged. For example, a reagent that will expire after midnight is loaded onto the instrument and scheduled to an assay before midnight of the expiration date. However, the test will not be completed until after midnight when that vial has expired. The instrument prevents any test runs with this vial in the future.

# Loading Plates

### Procedure

To load plates:

Step	Action
1	Insert each plate into its own transport frame with well A1 at the top left corner.
	<b>Note</b> : The transport frame has A1 imprinted onto the upper left corner to indicate the correct placement of a plate (as shown in the photograph below).
	<b>Note</b> : When using solid phase red cell adherence plates, verify that the plate is within its expiration date imprinted on the storage pouch and that the colored
	humidity indicator inside the storage pouch indicates that the plate has remained
	free from humidity.

Step	Action
2	Position the frame so that the guiding groove is on the right and the plate is inserted into the frame upright with well A1 at the top left corner. The barcode must be on the left edge of the plate.
	This image illustrates how the plate must be positioned in the transport frame.
	A BI
	Correct positioning of a plate in a transport frame
	A: Plate barcode
	B: Guiding groove
3	Make sure that the plate is correctly seated in the transport carrier depression and that it is not partially raised away from the bed of the transport frame.
	Warning: Incorrectly placing a plate in a transport frame may cause damage to the pipetting system and other modules on the instrument in addition to wasted resources.

Step	Action
4	Make sure that the required numbers of strips are inserted into the white plate frame, and that all strips are inserted in the correct orientation and pressed fully into position in the white frame.
	Warning: Loading a plate with an incorrect strip orientation results in invalid results and can create a biohazardous spill on the instrument. Incorrect orientation includes strips inserted upside-down in a white plate frame.
	<b>Warning</b> : Loading incorrect strips in a white plate frame can produce incorrect sample results without <b>Warning</b> . For example, loading Capture-R Ready-Screen strips in a barcoded Capture-R Select white frame is incorrect.
	Warning: When using Capture-R Ready-Screen or Ready-ID strips, load only the number of strips that are required for testing, as indicated in the <i>Resource</i> <i>Overview</i> window. Do not use Capture-R Ready-Screen or Ready-ID strips for testing that have previously traveled though the NEO Iris system on a plate frame, but were not used for testing.
	Note: The <i>Resource Overview</i> window displayed during the resource loading process confirms the number of strips needed, and the reader verifies the presence of the necessary number of strips before testing.
5	Open the <i>Plate Loading Tower</i> dialog by selecting the <b>Plates</b> button at the top of the <i>Resource Overview</i> window.

Step	Action
6	Open the door of the plate loading tower.
7	<text></text>
8	Ensure that the frame is inserted fully towards the back of the loading tower.
9	If there are no free slots, you may remove the plates in any slot with flashing green indicator lights. The light next to a free slot changes to continuous green to indicate that a new plate may be inserted at this position.

Step	Action
10	When you have loaded all required plates, close the loading tower door or press the <b>Scan new plates</b> button to initiate the plate barcode reading process.
	When the plate transport system is available, it takes each new plate from the loading tower and presents it to the barcode scanner for identification, then returns it to the same position in the tower.
	Warning: You must not place new plates into the tower while the transport system is accessing that module. Wait a few seconds after the transport has left the module, and then insert the new plates.
	Note: The plate barcode indicates the type of plate (such as blood grouping). This barcode determines the assays for which the plate may be used. The system displays these assays alongside the barcode information as shown below.
	Plate Loading Tower       No Status Plate ID       Assay Code       Stip Selection       Expiry Date       Process 4 >         15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       15       14       16       12       10       12       10 <td< th=""></td<>
	Attention: The plate tower door should be closed at all times, except when you are loading or unloading plates. This ensures that the plates are not accidentally disturbed.
	Note: The Cancel Schedule button can be used as a STAT resource interruption tool. If you have a STAT sample that you wish to begin immediately, but there are plates in the plate tower that are scanned and assigned into the system for other assays, but not yet started, you can remove all of the plates from the schedule using the Cancel Schedule button. This will allow more timely scheduling of the plate to be used for the STAT sample.
11	After loading a plate, make sure that the strip pattern in the <b>Plate Loading Tower</b> dialog ( <i>Strip Selection</i> tab) corresponds to the strips loaded in the plate.

#### Checking the Plate On Board Time

Although the NEO Iris provides no alerts if the on-board storage times for Capture strips are exceeded, there are tools available on the instrument to assist you in tracking this information.



**Warning**: You can load and store Capture strips (excluding Capture-R Ready-Screen and Ready-ID) on the instrument, outside of the storage pouch, according to the time **Limitations** printed in the relevant package inserts. However, it is the operator's responsibility to keep track of the on-board Capture strip storage time. The instrument will provide no alerts if the on-board storage time is exceeded and will not flag results generated from such strips.

To check the amount of time that the white plate frame has been on the instrument and when the white plate frame was first loaded:

Step	Action
1	Press the <b>Plate Loading</b> button on the Machine Monitor to access the Plate Loading dialog.
2	Select the plate ID that you are interested in by touching the Plate ID position on the <b>Plate ID</b> position on the dialog.
3	Select the <b>Strip Selection</b> tab. The Plate ID, Assay(s) to be run with this plate, the Use before date, the On Board Time, and the date and time that the plate was first loaded (First loading) are displayed at the bottom right of the screen.
	Plate Loading Tower         No Status         Plate ID         Assay Code         Stop Selection         Scan Plate         Assay Selection         Expiry Date         Processi (*)           Concel         15         14         12         3         4         5         7         8         9         10         11         12           Schedule         12         UA4310947         ABORH ABORH_E         10         0 </th

#### **Background Information**

- Strips that have been removed from the pouch and placed on the NEO Iris, but not used, should not be returned to the storage pouch.
- The On Board Time provided by the NEO Iris stops accumulating time once a plate is removed from the loading tower. When the plate is returned to the tower and the barcode is scanned, the time is not reset and the timer will resume counting. The time displayed represents the total "on-board" time.
- It is your responsibility to ensure that plates are not used beyond their approved on-board storage time. Each facility should implement procedures to ensure that this requirement is met.

### Starting Processing

#### When You Can Start Processing

You can start an assay run when all designated consumables for a given assay name on the *Resource Overview* window have green check marks.

#### Procedure

To start processing samples via the Resource Overview window:

Step	Action
1	Confirm that the consumable items designated for processing are selected and highlight the plate(s) that should be started.
2	plate(s) that should be started.         Press the Start button.         Stat         The software displays a running man animation on the screen while the tests are programmed into the schedule.         Image: State Start button.         Image: State
	well in this example. Under these circumstances, such wells will remain unused.
# Continuous Loading During Operation

#### Purpose

To maintain high throughput, you can access the instrument at any time to continuously load samples, replenish reagent, or empty waste while the instrument is processing.



**Warning**: Removing racks while the probes are accessing the tubes or vials from those rack results in damage to the probes and invalidated results. Wait until the LED for the rack is blinking green before removing.

### Loading Additional Samples or Requesting Additional Assays

You can load more samples if there are available loading positions, or you can request additional tests at any time on samples that are already loaded.

When all the samples in a rack have been pipetted and the LED for the rack is blinking green, you can remove the rack and insert a new one.

To load additional samples:

Step	Action
1	Open the <i>Start Run Assistant</i> by pressing the <b>Start</b> button on the Main Menu Bar.
	<u> </u>
	The system displays the Start Run Assistant dialog.
2	Click the Load Samples button to activate the Loading Bay dialog.

Step	Action		
3	You can remove a rack if the instrument has completed sampling the rack. These racks are marked by a flashing green LED in the loading bay lane.		
	<b>Warning</b> : The operator must first open the loading bay dialog in order to remove sample and reagent racks.		
	A continuous orange LED indic current processing. However, co outstanding tests remain for at as part of the current processir	ates that the samples are not required for the ontinuous orange can also indicate that : least one sample in the rack that is not scheduled ng.	
	Warning: Prematurely removing	g a sample or reagent rack will invalidate results.	
4	Load new samples and, if necessary, request that the assay be performed as described in <b>Loading the Samples</b> in this chapter.		
5	Click <b>Done</b> to close the <i>14-lane Bay</i> dialog	. The system redisplays the <i>Start Run Assistant</i> .	
6	If required, click the <b>Download Requests</b> button to obtain requested assays from the LIS.		
7	Click the Load Resources button. The system displays the Resource Overview window.		
8	Verify that you have loaded sufficient consumables onto the NEO Iris to complete the assays for the samples as described in <b>Completing the Sample Loading Process</b> in this chapter.		
	If T	Гhen	
	Sufficient consumables are available	Go to step 9.	
	Insufficient consumables are available	Replenish the necessary resource, and then go to step 9.	
9	Press the <b>Start</b> button.		
	Start The system generates a new schedule incom	rporating the new samples and assays.	

# **Reloading Plates**

After the system has completely used a plate, remove the plate from the instrument and insert a new one as needed.

To replace a plate:

Step	Action	
1	Remove the completed plates from the loading tower as indicated by the flashing green LEDs.	
2	Dispose of the used plates following standard laboratory practice.Marning: Used plates contain potentially biohazardous material. Wear protective gloves and clothing at all times when handling used plates. If any liquid is spilled, clean it up immediately following standard laboratory practice.	
3	Insert new plates into the Loading Tower as described in Loading Plates in this chapter.	

# **Replenishing Reagents and Controls**

Follow the procedure for Loading Reagents and Controls in this chapter.

The system displays a list of Missing (but necessary) reagents (refer to screen below).

	Revers III		🗰 🤣 🗐
Rack Area	Reagent IDs           Rack:         S           1         Reag1         203354102603599612           2         Image: Comparison of the second se	Reagent Propertie           Volume [ml]:         1           9.000         0           0.000         0           0.000         0           0.000         0           0.000         0           0.000         0           0.000         0           0.000         0           0.000         0           0.000         0	* nformation: ? ABORH, Anti-B, 2010-12-21 ??? ??? ??? ??? ??? ??? ??? ??? ???
Missing         Vol.▲           AnhA         0.1           AnhB         0.1           AnhD series 4         0.1           AnhD series 5         0.1           A1-Cell         0.1           One         Recall	5	10.000	<u>???</u>

## **Replenishing System Liquid**

To replenish the System Liquid, follow the procedure for **Filling the System Liquid** in **Chapter 10** - **Maintaining the NEO Iris**.

### **Removing Liquid Waste**

When the common waste container becomes full, the Waste icon in the Status Bar turns red and an audible alarm sounds.



A number of further pipetting and washer operations can occur (after the system displays the red waste icon) before the instrument stops pipetting and flags the results as invalid.

All liquid waste is collected in the common waste container. You can empty the liquid waste container at any time by draining the contents of the container into the shuttle container.



**Limitation**: If the common waste container is full and the waste icon turns red during sample processing, you will lose the interrupted testing run only if critical processing points exceed acceptable time limits because of excessive time removing the waste.

Examples of exceeding critical acceptable limits are:

- Too much time has elapsed between ABORh plate shaking and reading.
- A Capture-R Ready-Screen plate is incubated beyond the programmed time limit.

To drain the common waste container, follow the procedure for **Emptying the Common Waste Container** in **Chapter 10 - Maintaining the NEO Iris**.

# **Chapter 7: Test Results**

## In This Chapter

This chapter describes the methods for viewing, approving and exporting test results.

CHAPTER 7: TEST RESULTS	7-1
Accessing the Results Screen	7-2
Sample View and Plate View Icons and Symbols	7-7
Using Tool Tips	7-9
Viewing Test Details	7-10
Approving Test Results	7-17
Exporting Test Results	7-18
Viewing Archives	7-20

# Accessing the Results Screen

The Results screen is where test results can be viewed.

To access the Results screen:

Step	Action
1	Click the <b>Results</b> button on the Main Menu Bar of the Machine Monitor.
	Results
2	The <i>Results</i> screen with the most recently used screen view (sample or microplate) appears.
	Image: Sample :       Assey :       All result       Mapproval Status       Actions         View       Plate :       All result       Assey :       All result       Plate :       Assey :       Assey :       All result       Plate :       Plate :       Assey :       Assey :       All result       Plate :       Plate :       Assey :       Assey :       Assey :       Assey :       All result       Plate :       Plate :       Plate :       Assey :       Assey :       Assey :       Assey :       Assey :       Plate :       Assey :       Assey :       Plate :       Assey :       Assey :       Plate :       Assey :       Plate :
	2009/08/07 13:44:42       67/78       Syph       Pending         2009/08/07 13:44:42       P016002P       Syph       Pending         2009/08/07 13:44:42       ME5987515       Syph       Pending         2009/08/07 13:44:42       ME5987571       Syph       Pending         2009/08/07 13:44:42       ME5987571       Syph       Pending         2009/08/07 13:44:42       ME5987571       SYPH_RPT       Pending         No Sample Results matching search criteria       ME5987571       SYPH_RPT       Pending         No Sample Results matching search criteria       ME5987571       SYPH_RPT       Pending         Database Name: Galileo-PC/Galileo_SOL       OK

# Viewing Results

To view results:

Step	Action	
1	Click the <b>Sample View</b> button in the View area of the window.	
	- Sample View (see step 2)	
	- Plate View (see step 3)	
	You can only select one of these buttons at a time. The selected button exhibits a "pressed"	
	appearance. The Results area defaults to the most recently used of these two views when first	
	selected.	
2	The Sample View Results screen appears.         Very       Deta Selction         Plate:       All pets:         All pets:       Plate:         All pets:       Plate:         Donation:       All pets:         Donation:       All pets:         Donation:       All pets:         Plate:       Plate: <th< th=""></th<>	
	Database Name: Galileo-FC/Galileo_SQL	
	The Data Selection controls affect which records are displayed in the list.	

#### Chapter 7: Test Results

Step	Action
3	The <i>Plate View Results</i> screen appears.
	View       Data Selection         Plate :       All         Sample :       All         Donation :       All         To :       2009/09/17
	Clate         Plate Read Date         Plate ID         Assay Name         Tests         Flags         Edited         Approved         Exported         Comment           2009/09/17 10:06:55         2009/09/17 10:06:35         UA5364570         mAgSor37         91         Plate has interpreted without           2009/09/17 09:19:19         2009/09/17 08:10:39         2009/09/17 08:10:39         2009/09/17 08:10:24         VA5364580         46           2009/09/17 08:10:39         2009/09/17 08:10:24         UA5364580         QCTEST         8         Plate has interpreted without           2009/09/17 08:10:39         2009/09/15 08:50:25         UA5415035         ABORH         2         Plate has interpreted without
	2009/09/15 14:26:03         2009/09/15 14:26:03         DP03501007         ExtendOP         1         Image: Comparison of the Comparison of t
	2009/09/14 14:00:43         2009/09/14 14:00:26         UA5415034         FWD_ABORH         4         Plate has interpreted without           2009/09/14 13:54:55         2009/09/14 13:54:55         UA5415035         FWD_ABORH         12         Plate has interpreted without           2009/09/14 13:54:55         2009/09/14 13:54:57         UA5415085         FWD_ABORH         12         Plate has interpreted without           2009/09/14 13:34:55         2009/09/14 13:34:55         UA5415085         FWD_ABORH         12         Plate has interpreted without           2009/09/14 13:34:55         2009/09/14 13:34:56         UA5415085         FWD_ABORH         12         Plate has interpreted without           2009/09/14 13:34:51         2009/09/14 13:34:51         UA5415085         FWD_ABORH         12         Plate has interpreted without           2009/09/14 13:34:59         2009/09/14 13:34:51         UA5415085         FWD_ABORH         12         Plate has interpreted without           2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/14 13:34:59         2009/09/
	Database Name: Galileo-PC/Galileo_SQL
	The Data Selection controls affect which records are displayed in the list.

Step	Action	
4	For both the Sample and Plate Views, enter any of the following criteria:	
	• In the Plate field, type a plate ID. This field defaults to <b>All</b> .	
	• In the Sample field, type a Sample ID. This field defaults to <b>All</b> .	
	• In the Donations field, type a Donation ID. This field defaults to All.	
	<ul> <li>From the Assay drop-down list, select the assay from the list of available assays. You can select All from the list to include all assays.</li> </ul>	
	<ul> <li>Select the From and To date drop-down lists to select a different from or to date.</li> <li>Select a date from the displayed calendar box.</li> </ul>	
	• From the Approval Status area, select an approval status. This area defaults to <b>All Tests</b> . Other options are <b>Not Reviewed</b> , <b>Approved Only</b> and <b>Not Approved</b> .	
	• From the Result Types area, select a result type. This area defaults to <b>All Types</b> . Other options are <b>Unusual</b> and <b>Normal</b> .	
	If you change any of the data selection controls, the list updates immediately.	
	The <b>Other Filters</b> button is used to access the Other Filters dialog. This dialog allows additional criteria for filtering results listed on the results screen. Two (2) tabs are available,	
	<b>Reagents</b> and <b>Miscellaneous</b> . The Reagents tab allows filtering via <b>Reagent Name</b> and/or	
	<b>Lot Number</b> . A drop down list is used for both, where selections can be made. The Miscellaneous tab allows filtering via <b>Instrument ID</b> . A drop down list is used to make a	
	selection. Selecting the Miscellaneous tab allows an instrument identification to be selected	
	from a list of identifications saved in the database. The default selection for all of the drop	
	down boxes is <b>All</b> .	
	Instrument ID:         Image: Construction of the second seco	
	Use the Cancel button of the Other Filters dialog to cancel the selections made or use the	
	<b>OK</b> button to apply the filtering selections made. If you change the other filters criteria and apply them using the <b>OK</b> button, then the result list updates immediately.	

Step	Action	
	Other Filters     X       Reagents     Miscellaneous       Reagent Name :     All       Lot Number :     Instrument ID :	
	OK Cancel OK Cancel	
	Data identifier columns on both the plate view and sample main view results display have adjustable widths.	

# Sample View and Plate View Icons and Symbols

The table below contains descriptions of the icons/symbols used in the *Sample View* and *Plate View* windows:

Icon/Symbol	Description
Details	The <b>Details</b> button displays the <i>Test Details</i> dialog which contains information for the item that was selected in the main <i>Results</i> screen.
Approve	The <b>Approve</b> button displays the <i>Approval</i> dialog which contains approval information for the item that was selected in the main <i>Results</i> screen.
Print	The <b>Print</b> button displays the print dialog. This enables you to print results. Refer to <b>Chapter 8 – NEO Iris Reports</b> for details on how to print reports.
Export	The <b>Export</b> button displays a dialog that has a list box containing the currently approved, but not exported, tests. This enables you to export results if a functioning interface is set up.
	Refer to <b>Attachment 2</b> for the interface configuration requirements.
Edit	The <b>Edit</b> button displays the <i>Test Details</i> dialog. Refer to <b>Chapter 8 – NEO Iris Reports</b> for a list of assays with equivocal reactions that cannot be edited.
X Void	The <b>Void</b> button displays a confirmation dialog that displays a message signaling that continuing will void a chosen test result.
	Pressing the <b>Delete</b> key on the keyboard has the same effect as the <b>Void</b> button. A voided test will not appear in later search results.
<mark>₽</mark> efresh	The <b>Refresh</b> button allows you to update the list of results on the screen while the instrument is running, and without the need to exit and then reenter the <i>Results</i> screen. This button becomes active when new results become available.
Δ	In the Flags column, the system displays the warning symbol if the sample has any associated warnings. You must review any test results that have an associated warning flag.

Icon/Symbol	Description
0	The system displays the error symbol if the sample has any associated errors. The system also displays this icon to indicate non-critical process incidents that occurred during sample processing that do not affect the sample result.
	In the Edited column, the system displays the edited symbol if any results for that test have been edited.
۸	The software has built-in parameters to automatically designate a result for further testing. This function is called Validate Method. It is represented by the <b>Validate</b> icon. An example of this designation is qualifying all initially Rh (D) negative tests as Weak D Pending. Weak D testing subsequently resolves the sample result. These parameters are configurable on an instrument-by-instrument basis.
<b>√ ×</b>	In the Approved column, the system displays the green check mark symbol if the test has been approved. The system displays the red <b>X</b> symbol if the test has been unapproved. The system displays no symbol if the test has not yet been approved. In the Exported column, the system displays the green check mark symbol if the test has been exported. The system displays no symbol if the test has not yet been test has been exported.

# Using Tool Tips

When you leave the mouse cursor stationary over a line, error symbol, or warning symbol for one-half second or more, the system displays a tool tip. If the cursor is held over the sample ID, then the patient information tool tip appears. If the cursor is held over an error/warning symbol, then the ID and error information appear.

The screen image below displays an example of the tool tip.



A flashing screen can be displayed if the cursor is held over a symbol that has an excessive amount data linked with it. The tool tip associated with that particular symbol cannot be displayed because the size of the tool tip, determined by the amount of data linked with it, exceeds the size of the screen display.

# Viewing Test Details

### Procedure

To view test details:

Step	Action
1	Click the <b>Details</b> button in the <i>Actions</i> area when in the <i>Sample View</i> mode.
2	The system displays the <i>Test Details</i> dialog, which contains the <i>Test Overview</i> tab for the first test that was selected when you clicked the <b>Details</b> button. The completion of the <i>Test Overview</i> tab fields is optional. Alternatively, double-clicking a line in the list window box will also display the <i>Test Details</i> dialog.
	<ul> <li>This screen consists of 4 additional tabs:</li> <li>Results</li> <li>Reagents</li> <li>Event Log</li> <li>Plate Views</li> <li>Mote: If the Plate View mode is used, then the default tab that the <i>Test Details</i> dialog displays is <i>Plate Views</i>.</li> </ul>

### **Results** Tab

The left side of the dialog, below the *Results* tab, is a list box that displays all samples for the selected plate. From this list box you can select the particular sample that you want to view.

The right side of the dialog displays test details for the selected sample. On the far right-hand side of the dialog, you can make edits to the results.

To edit a well reaction, choose one of the available results from the *Revised Grade* drop-down menu.



<u>Warning</u>: Refer to **Attachment 1 for NEO Iris Operator Manual** for a list of assays with equivocal reactions that cannot be edited.

Ы	est Details - Plate: R	73403539, Sampl	e: W053316	i821148, Assay: 3_	Cell						x
ſ	Test Overview Resul	ts Reagents Eve	ent Log   Plai	te Views							
	Sample No.	Edited	Well			Well	Reagent	Well	Reaction	Original Grade	Revised
	W053316821148	Negative	A02D02		<b>^</b>	100	c-11.4	Indge	2 Cal		
	W053316001677	Negative	E02H02			AU2	Cell I		<b>O</b>	$\bigcirc$	
	W053316604447	Negative	A03D03			B02	Cell 2		(3)	(-)	- 🔽
	W053316050318	Negative	LU3HU3 A04 D04				- 1 -		Ă	X	
		Negacite	Homeon			CU2	Cell 3		9	$\bigcirc$	
						D02	IgG Pos		(83)	(3)	3 🔻
								L Traciti	$\sim$	$\sim$	
					-						
							Negative		] _⊺e	st Status	
							Inegative		1	Approv	red
									1	Report	ed
					Co	omment :				Export	ed
						1			<b>4</b>		
						<b>V</b>			Drint	-	
1					_	whhlone	Export		Princ		liose

#### **Reagents Tab**

The *Reagents* tab displays details for all reagents that were used for the currently selected results.

Test D	etails - Plate: l	JA5682662,	Sample: I	LQ136083, <i>i</i>	Assay: ABORH				×
Test (	Overview Resul	Its Reagents	Event Lo	og   Plate View	vs				Φ
Well	Reagent	Lot Number	Qualified	Expiration	Last QC	QC Expires	Barcode	Edited	Edited By
A02	Anti-A	101711	1	2011/03/23	2010/02/25 08:01:47	2010/02/26 08:01:32	101082117110949625		
B02	Anti-B	203281	1	2011/05/07	2010/02/25 08:01:47	2010/02/26 08:01:32	203127112819824198		
C02	Anti-D series 4	504732	~	2011/03/07	2010/02/25 08:01:47	2010/02/26 08:01:32	504066117321675430		
D02	Anti-D series 5	505571	1	2011/05/07	2010/02/25 08:01:47	2010/02/26 08:01:32	505127115711037806		
E02	Mono ctrl	492051	$\checkmark$	2011/08/06	2010/02/25 08:01:47	2010/02/26 08:01:32	492218110511853900		
F02	A1-Cell	111737	1	2010/02/26	2010/02/25 08:01:47	2010/02/26 08:01:32	111057107370326511		
G02	B-Cell	113737	1	2010/02/26	2010/02/25 08:01:47	2010/02/26 08:01:32	113057107370185817		
•									•
							3		
							Print	C	lose

#### **Event Log Tab**

The *Event Log* tab displays the detailed list of events that occurred for testing of a given sample or plate. If the plate has been used for more than one run, the event log lists the events for the run currently being viewed.

	Trace Data														1	
	2010/02/25	12:28:15	"LQ13	5083 foun	d at track	1, pos 1 in	Sample B	ay"								
	2010/02/25	12:28:15	"10108	"101082117110949625 found in 14-lane Bay at track 14, pos 2"												
	2010/02/25	12:28:15	"20312	"203127112819824198 found in 14-lane Bay at track 14, pos 1"												
	2010/02/25	12:28:15	150406 UE0E12	"504066117321675430 found in 14-lane Bay at track 14, pos 3"												
	2010/02/25	12.20.10	1110P	"b0b12/115/1103/806 found in 14-lane Bay at track 14, pos 4"												
5	2010/02/25	12:28:15	111305	"11105/10/3/0326511 found in 14-lane Bay at track 14, pos 6" "1120672107307195917 found in 14-lane Bay at track 14, pos 6"												
j I	2010/02/25	12:28:15	"49221	"492218110511853900 found in 14-lane Baulat track 14, pos 7												
2 2010/02/25 12:26 To San Traditional Traditional Bay But tack 14, pp 5																
2	2010/02/25	12:28:42	28.42 plate moved to Reader													
	2010/02/25	12:28:53	Camera Reader results stored in file `D:_4b86b355_0.bbx`													
į.	2010/02/25	12:28:53	StripCle	ean Read	(Reader)-	Results:										
	2010/02/25	12:28:53	546	857	850	845	855	848	851	859	858	860	854	871		
	2010/02/25	12:28:53	467	845	840	841	848	848	847	853	850	852	863	869		
	2010/02/25	12:28:53	514	843	841	841	844	839	843	843	848	847	852	857		
1	2010/02/25	12:28:53	468	844	850	847	841	844	846	852	847	847	856	862		
	2010/02/25	12:28:53	524	840	839	834	842	846	839	849	847	851	856	872		
	2010/02/25	12:28:53	536	836	839	841	841	842	839	843	835	841	842	860	1	
UDSUT	2010/02/25	12:28:53	536	836	839	841	841	842	839	843	835	841	842	1	860	

When viewing a plate that has failed controls and has not been saved, it is still possible to view the **Event Log** and **Plate Views** tabs. However, in this instance, the **Test Overview**, **Results** and **Reagents** tabs cannot be accessed and the tabs will be grayed out, as shown in the image below.

Test	Details - Plate: SC132	01332, Assay: WEAK_D	x
Tes	t Overview Results Re	agents Event Log Plate Views	
	Trace Data		
ess Record Results Log	Trace Data 2010/02/03 14:22:17 2010/02/03 14:22:17	"R142959 found at track 4, pos 1 in Sample Bay" "R142960 found at track 4, pos 2 in Sample Bay" "R142962 found at track 4, pos 3 in Sample Bay" "R142965 found at track 4, pos 5 in Sample Bay" "R142965 found at track 4, pos 5 in Sample Bay" "R142967 found at track 4, pos 7 in Sample Bay" "R142967 found at track 4, pos 7 in Sample Bay" "R14297 found at track 4, pos 7 in Sample Bay" "R14297 found at track 4, pos 7 in Sample Bay" "R14297 found at track 4, pos 8 in Sample Bay" "R14297 found at track 4, pos 9 in Sample Bay" "R14297 found at track 4, pos 9 in Sample Bay" "R14297 found at track 4, pos 10 in Sample Bay" "R14297 found at track 4, pos 11 in Sample Bay" "R14297 found at track 4, pos 11 in Sample Bay" "R14297 found at track 4, pos 12 in Sample Bay" "R142987 found at track 4, pos 13 in Sample Bay" "R142981 found at track 4, pos 14 in Sample Bay"	
og Proc	2010/02/03 14:22:17 2010/02/03 14:22:17 2010/02/03 14:22:17 2010/02/03 14:22:17 2010/02/03 14:22:17	"R142984 found at track 4, pos 16 in Sample Bay" "R142987 found at track 5, pos 1 in Sample Bay" "R142990 found at track 5, pos 2 in Sample Bay" "R142991 found at track 5, pos 2 in Sample Bay"	
Instrument L	2010/02/03 14:22:17 2010/02/03 14:22:17 2010/02/03 14:22:17 2010/02/03 14:22:17 2010/02/03 14:22:17	"R142391 found at track 5, pos 4 in Sample Bay" "R142992 found at track 5, pos 4 in Sample Bay" "R142998 found at track 5, pos 5 in Sample Bay" "R142998 found at track 5, pos 7 in Sample Bay"	-

# **Plate Views Tab**

The *Plate Views* tab displays several options to view plate information.

iam	ple No.	Edited ABO	Rh	Well
	R132522	В	Positive	A02H02
	R132493	0	Positive	A03H03
	R132559	0	Positive	A04H04
	R132521	А	Positive	A05H05
1	R132488	0	Positive	A06H06
	R132554	0	Positive	A07H07
	R132547	А	Positive	A08H08
	R132532	А	Positive	A09H09
	R132517	в	Positive	A10H10
	R132494	A	Positive	A11H11
	R132527	0	Negative	A12H12

To view Plate Views from this tab:

Step	Action
1	Click the plate view +/- button.
	The system expresses the results in plate wells as + or
	$\odot \bigcirc$

Step	Action
2	Click the <b>Grade</b> view button to display the same plate with results expressed as a graded numeric value for each reaction.
	Grade
3	Click the <b>Reaction Strength</b> view button to display the same plate with results expressed
	numerically.
	94 Reaction Strength
	Any value over 99 will be displayed as 99. All numbers are rounded so it is possible to see
	the same value classified as positive and negative.
4	Click the <b>Image</b> view button to display the same plate with the well images shown

Manually edited wells will be flagged with an edited symbol **II**.

Some icons/symbols are common to most Plate Views, with the single exception of the Image icon.

The table below contains descriptions of these common symbols.

Symbol	Description
<b>(</b> )	The Equivocal symbol indicates an equivocal result for a given well.
×	The Invalid Well symbol indicates an error status for a given well, such as quantity not sufficient (QNS) or other error.
0	The Blank symbol indicates wells that are not occupied.

#### Viewing Additional Sample Results

If a sample entry is highlighted in the Sample View window and a previous result exists for that same sample with the same assay (such as a previous blood group NTD now resulting with a defined blood group), and you press the **Details** button, the system displays a dialog warning you that other results

exist. The *Additional Sample Result(s) located* dialog offers you the option to proceed if you press the **OK** button. Press the **Cancel** button to cancel.

Additio	onal Sample	Result(s) loc	ated	
<b>(1)</b>	Another Resu	It exists for thi	s Sample for the sa	me Assay.
Ŷ	Select Cancel Select OK to (	to view Result continue and vi	(s) by re-selecting I ew Results related	from the Sample list. to this plate.
		ОК	Cancel	l

# **Exiting Test Details**

When you are finished, click the **Close** button of the *Test Details* dialog.



The system returns to the *Results* screen.

# Approving Test Results

To approve test results:

Step	Action
1	Select the plate or sample results that you would like to approve in the main <i>Results</i> screen. Click the <b>Approve</b> button in the <i>Actions</i> area of the screen.
	$\underbrace{\text{Note}}_{\text{a result.}}$ The <i>Approve</i> dialog can also be selected while in the <i>Test Details</i> dialog of a result.
2	The system displays the Approve dialog.         Image: No.         Image:
3	Click the <b>Approve</b> button in the <i>Approve</i> dialog to approve all tests where check boxes are checked.

# Exporting Test Results

To export test results:

Step	Action							
1	The system displays a <i>Result Approval</i> dialog after a result is approved. This dialog gives you the option to export the approved results to the LIS.							
	Would you like to export the approved tests now?         Yes       No         Cancel							
	Press the <b>Cancel</b> button of the <i>Result Approval</i> dialog to cancel the approval process.							
	Note: By pressing the <b>No</b> button, the results are approved but not exported, and are then subsequently available for export, if required, from the main results screen. In this instance, the approved results can be exported by selecting the <b>Export</b> button from the <i>Action</i> s area of the main results screen when those results are selected.							
	Export							
	Before approval, the <b>Export</b> button is not available for use when unapproved results are							
	selected in the main results screen.							
	Export							

Step	Action
2	Press the <b>Yes</b> button of the <i>Result Approval</i> dialog to export. The system displays the <i>Exporting</i> dialog with the correct report configured and the selected samples or plates listed.
	Exporting Search Type: Report Name : Galileo Export Sample - Current Results
	Sample No.         Donation No.         Assay         Result 1         Result 2         Flags         Edit           N48219         5055         IgG_XM         IgG Comp (Check ABO Comp)         IgG
	✓ Manakati Manaka
3	Select the individual results to be exported using the line check boxes or select all results using the <b>Select All</b> button of the <i>Exporting</i> dialog. Press the <b>Cancel</b> button to cancel the export process if required.
	Note: Verify that all samples to be exported have been selected on the <i>Exporting</i> dialog after using the line check boxes or selecting all results by pressing the <i>Select All</i> button.
4	Press the <b>Export</b> button of the <i>Exporting</i> dialog to export. The monitor display returns to the previous screen that was used before the <i>Exporting</i> dialog, following the result export.

# **Exiting Test Results**

To exit *Test Results*, press the **OK** button. The monitor display will return to the Machine Monitor screen.

# Viewing Archives

This section describes how to view archives.

### Procedure

To view archives:

Step	Action
1	If result archives are stored on disks, insert a previously archived DVD+R disk into the disk drive.
2	Select the <b>Login</b> button on the main menu bar to display the <b>Login</b> dialog.
3	Enter your Name and Password into the Login dialog. Select the checkbox for Use DMS Archive and then press the OK button of the Login dialog. NEO Iris - Login Please enter your name and password Name : Password : Start Action Restore Databases from Backup Use DMS Archive OK Shutdown The list of all archives performed on the instrument will then be displayed on the Select DMS Archive dialog. Note: The Select DMS Archive dialog will display all archives ever performed on the particular instrument. The most recent archive will be at the bottom of the list.

Step	Action
4	Select the archive you would like to view from the displayed list of all archives performed on the instrument. The selected archive must match the archive stored on the archive medium attached to the PC (or DVD+R disk inserted into the drive in step 1 above).
	Name         Comments         DB Type           20090709_20090806-40         SQL Server           20090709_20090806-5         SQL Server           20090709_20090806-6         SQL Server           20090709_20090806-7         SQL Server           20090709_20090806-8         SQL Server           20090709_20090806-9         SQL Server           20090709_20090806-9         SQL Server           20090709_20090806-9         SQL Server           20090709_20090818         SQL Server           20090818_20090818         SQL Server           20090819_20090819         SQL Server           20090819_20090901         SQL Server           20090901_20090901         SQL Server
	20090911_20090911       SQL Server         20090929_20090930       SQL Server         20090929_20100120       SQL Server         20090929_20100126       SQL Server         Image: State OK button and allow the software to retrieve the archive files selected and do not
	Press the <b>OK</b> button and allow the software to retrieve the archive mes selected and do not press any keys during this process. <b>Note</b> : When the archive is viewed, the software will first create a copy of the archive on the hard drive. The local copy will be opened for viewing; the archive on the disk is never opened. This reduces the potential for corrupting the archive.
5	Select the <b>Results</b> button on the Main Menu Bar to display the archived results on the results screen.
	<b>Note</b> : You can customize your viewed result data, for example by date range, by using the methods described in this chapter.
6	Press the <b>OK</b> button of the results screen to close the screen after you have finished viewing the archived results.
7	To target the instrument software back to the list of current results, select the <b>Login</b> button on the Main Menu Bar to display the <b>Login</b> dialog.

Step	Action
8	Enter your <b>Name</b> and <b>Password</b> into the <b>Login</b> dialog. Do not select the checkbox for <b>Use</b> <b>DMS Archive</b> and then press the <b>OK</b> button of the <b>Login</b> dialog to open the working database.
9	If the archive was viewed from a disk, remove it from the PC disk drive and store it in a safe place, according to the recommended storage conditions of the disk manufacturer.

# **Chapter 8: NEO Iris Reports**

## In This Chapter

This chapter describes the NEO Iris reports.

CHAPTER 8: NEO IRIS REPORTS	8-1
Reports Overview	8-2
Parts of the Report	8-4
Plate Based Reports	8-10
Sample Based Reports	8-13
Current Reports	8-15
Quality Control Reports	8-20
Reagent Reports	8-23
Accessing Plate Based Reports	8-25
Accessing Sample Based Reports	8-27
Printing Reports	8-29
Printing Reagents Report from Test Details	8-35
Test Results and Interpretation	8-37

# Reports Overview

#### **Result Report Types**

NEO Iris reports follow a general formatting scheme that is described by the following three main types. Each result report is assay specific and the different assay reports vary in the format details.

So that you can obtain all necessary printed information, the three main types of result report are:

Type of Report	Description
Plate Based	This type of report is plate based, and can only report sample results generated on a single plate run. This report type displays samples in the order they appear on the plate and includes control information. Plate based reports are the only method to print control results.
Sample Based	This type of report type displays the sample(s) selected. Samples from different plate runs can be printed on one page as a batch report. Samples are printed in alphanumeric order according to their plate ID. Only sample interpretations are printed on this report; there are no well results printed. If you select samples from an assay that is not the same as the selected report assay type, the report displays blank lines where these samples would have been reported.
Current Result	This type of report displays only the most current results for a sample. This type of report includes red blood cell antibody identification assay reports. Only one sample is listed per page along with the reagents used (not including the control information). This type of report includes Group and Screen reports.



**Note**: The report templates used in this chapter are examples chosen to illustrate the various report types and their respective design. Please, contact your local Immucor Representative to get more information about the reports that are available.

#### **Generating Reports**

A maximum of two reports can be configured to automatically output when the system interprets the result file. Immucor configures this for printed reports or an export file to the LIS. The export file to the LIS is defined as a report output.

You can also manually print reports in the results section of the software.

# Parts of the Report

# **Example of a Report**

This example of a report exhibits features common to all reports and is described area-by-area in this chapter.

			3_	Jell F	vebor	L			IRIS	
Facility Name: Address:	Immucor R&D Bidg 7000 Norcross, GA		Assay ID Operator Date Re: Time Re	: 3, ID: In Id: 21 Id: 11	_Cell ndsayw 016/06/21 8:25:32 P	Ins So Aa M Re	strument SN: Aware Version: say Version: sport Version:	5030090 1.6.11.1 4.00 05 3.00	1873 2. Nov12 3.	11 SP1 00 30Apr
Plate ID: Lot Number: Expiry Date:	R73400111 R734 2016/07/15									
Sample ID: LQ258281		Interpretation: Negative	Cell	Cell I	Cell III	Pos Ctrl 4	Current Interp Negative	retation:	Reference R73400111	Flag
				-						

The report consists of three (3) sections:

- ► Report Header
- ► Report Body
- ► Report Footer

# **Report Header**

The following is an example of a report header.

імми	JCOR	3_Cel	ll Report		NEO IRIS	+
Facility Name: Address:	Immucor R&D Bldg 7000 Norcross, GA	Assay ID: Operator ID: Date Read: Time Read:	3_Cell lindsayw 2016/06/21 18:25:32 PM	Instrument SN: Software Version: Assay Version: Report Version:	5030090873 1.6.11.1 4.00 05Nov12 3.00	2.11 SP1 3.00 30Apr14
Plate ID: Lot Number: Expiry Date:	R73400111 R734 2016/07/15					

The report header consists of the following:

Graphic	Name	Description
IMMUCOR.	Immucor Logo	The Immucor logo is located in the top left corner of all NEO Iris reports.
Facility Name: Immucor R&D Address: Bldg 7000 Norcross, GA	Facility	The blood testing facility from which you generate reports and its address are located immediately below the Immucor logo.
	Report Title	The report title is located at the top center of all reports. This title identifies the type of report and what it represents.
Assay ID: 3_Cell Operator ID: lindsayw Date Read: 2016/06/21 Time Read: 18:25:32 PM	Details	The specific run details for the reported results are located under the report name, which includes the instrument plate read date and time. This corresponds to the information in the <i>Plate Read Date</i> column of the results display. <u>Note</u> : Group and screen reports do not list any information in this area. Reagent reports do not list the Operator ID.
NEO + IRIS	NEO Iris Logo	The NEO Iris logo is located in the top right corner of all reports.

Instrument SN: 5030090873 Software Version: 1.6.11.1 Assay Version: 4.00 05Nov12	Instrument Serial Number, Software Version, Assay Version	The instrument serial number, software versions, and assay versions related to the reported results are located under the NEO Iris logo. If the report lists results for multiple assays, the report identifies assay versions for each assay. If it is a Reagent Report, no assay version is included.
Report Version: 3.00	Report Version Number	The report version number is located directly below the software and assay version numbers, at the right hand side of the report header.

# **Report Body**

Plate ID:         R734001           .ot Number:         R734           Expiry Date:         2016/07/	111							
Sample ID:	Interpretation:	Cell I	Cell II	Cell III	Pos Ctrl	Current Interpretation:	Reference	Flag
LQ258281	Negative	-	-	-	4	Negative	R73400111	

The report body contains the reported results and Plate details.

Each report includes a processing **Flag** column within the report result grid. This column is either left blank or is populated with a 2 or a 3. The key states that 2 indicates a **Warning** and 3 indicates an **Error**. A blank table cell indicates no associated flags. Examples are shown below.

erence	Flag	nce	Flag
5114059	2	3156	3
5114059	2	3156	3
5114059	2	3156	3
5114059	2	3156	3



**Note**: **Warning** (2) flags do not require you to review any associated test results. There are three scenarios that generate a 2 flag. These scenarios are described in the table below:

Warning Text	Interpretation	Potential Issue
ID entered manually	A reagent, sample or plate ID was entered manually.	Loading bay or plate barcode scanner cannot read a barcode. Reagent, sample, or plate barcodes cannot be read due to poor quality of barcode.
No barcode read	The plate barcode scanner was unable to read the barcode on the plate during assay processing.	Plate barcode scanner does not read.
Warning in device	This is the warning assigned to any result where plate jam errors were recovered from successfully. If error recovery was unsuccessful or untimely, then associated results will be invalidated.	Error recovery after plate jam.

**Error** (3) flags are either associated with results where the interpretation is invalid or with unused wells. Neither situation requires you to review any associated test results.

Other report features are discussed later in this chapter.

#### **Report Footer**

The following is an example of the report footer.

Key:	*INV* = Invalid Interpretation No_Int = No Interpretation Ctrl Fail = PC < 3+ or NC > 0	- 1-4 ? X	= Negative Result = Graded Positive Result = Equivocal Result = Invalid Result	Flags:	2 3	= Warning = Error		
2016/0	6/21	18	8:25:59 PM				Page:	1

Reagent information is listed in a separate grid for some reports.
Accov Nama:	Reagent Name:	Lot Number:	Lot Expiry:
Assay Name.	Reagent Name.	Lot Number.	Lot Expiry.
ABORH	Mono ctrl (10ml)	492241	2018/03/01
ABORH	A1-Cell (10ml)	111083	2016/07/08
ABORH	Anti-A (10ml)	101746	2017/09/01
ABORH	B-Cell (10ml)	113083	2016/07/08
ABORH	Anti-B (10ml)	203620	2017/10/20
ABORH	Anti-D series 4 (10ml)	504715	2018/03/03
ABORH	Anti-D series 5 (10ml)	505741	2017/08/11
3_Cell	LISS (10ml)	211750	2017/03/29
3_Cell	Cap-R Ind Cells (10ml)	221650	2016/07/07

At the bottom of most reports is a key with symbols to assist you in interpreting the reports. An example is shown below. The key is populated with information pertinent to the report upon which it is printed. Therefore different types of reports will be printed with different keys. For example, the crossmatch assay includes **Comp = Compatibility** in the key. This item in the key is only found for crossmatch assay. Quality control reports do not print with a key included.

Key:	*INV* = Invalid Interpretation	1-4= Graded Positive Result	Flags:	2 = Warning
	No_Int = No Interpretation	? = Equivocal Result		3 = Error
	NTD = No Type Determined	X = Invalid Result		
	Ctrl Fail= PC < 3+ or NC > 0	= Negative Result		

The date and time that the report is generated, as well as the page number, are also located at the bottom of the page. An example is shown below.

2016/06/22 16:11:01 PM Page: 2

## Plate Based Reports

#### **Report Details**

The details of the plate used to determine the reported results are located on the left side of the report, above the results grid. The results grid is in the center of the report. The column titles vary slightly, depending on the assay reported.

Facility Name: 1 Address: 1	mmucor R&D Bidg 7000 Norcross, GA R75400111	Assay ID Operator Date Res Time Res	: 3, ID: lin id: 20 id: 18	,Cell dsayw 016/06/21 0:25:32 Pf	in: So Aa M Re	strument SN: 503 Aware Version: 1.6 say Version: 4.0 port Version: 3.0	00090873 111.1 0 05Nov12 0	2.11 SP1 2 3.00 30Ap
Lot Number: Expiry Date: 2	R734 2016/07/15							
Sample ID: LQ258281	Interpret	ation: Cell	Cell	Cell III	Pos Ctrl	Current Interpreta	tion: Refe	erence Flag
C. C				-		Negative	PLTD	400111
			_				_	
			-					
			_					
			_					
							_	
Key: "INV" No_Int Cbi Fail	= Invalid Interpretation = No Interpretation = PC < 3+ or NC > 0	- = Ne 1-4 = Gri ? = Eq	gative Re aded Pos uivocal R	esult itivo Result lesult	ult	Flags: 2 3	= Warning = Error	2
		X = inv	and Hes	UFL.				

**Note**: The printed plate based report option lists the sample IDs in the pipetted plate order. It is important to note that the pipetted plate order is not the same as the sample order in the 14-lane Bay or the order the samples were pipetted to the plate. This pipetting technique optimizes the pipetting procedure and overall workflow. The controls are listed on the plate reports. When printing from plate view, only plate-based reports are available.

## **Report Columns**

The table below describes the Plate Based Report columns.

Column	Description
Sample ID	The patient identification barcode information as it appears on the tube.
Interpretation	The instrument interpretation from the original results created on this plate run.
Result text	The result text column(s) represent each test well (with a description of the well, such as Cell I or Anti-A) and the result that the system used to determine the interpretation.
Current Interpretation	This interpretation may be the same as the original interpretation. If you repeated the sample with the same assay and had a different result, edited the sample, or the sample was run with the validate method, this interpretation is different than the original interpretation.
Reference	This column indicates where the current interpretation is from. If you have not adjusted the interpretation from further testing (such as VALIDATE), this ID is the original plate ID. Otherwise, it represents the plate that the current interpretation was generated from (such as a second run of the same sample with the same assay).
	If the current interpretation is from an edit, the report displays 'Manual Edit' in this column. If the current interpretation is from the result of a validate method, the report displays 'VALIDATE' in this column.
Flag	This processing flag column is either left blank or is populated with a 2 or a 3. The key states that 2 indicates a <b>Warning</b> and 3 indicates an <b>Error</b> . A blank table cell indicates no associated flags.

## Sample Based Reports

#### Sample Based Report Types

Sample based reports include reports for samples tested by a single assay.

#### Sample Based Report Details

The sample based report lists all highlighted samples for one assay. This repeat can print samples run on different plates and will report the samples in alphanumeric order. Refer to **Reagent Reports** later in this chapter for information regarding reagent reports.

Sample ID:         Inferpretation:         Reference         Plate ID:         Date Rest:         Time Read:         Operator ID:         Flag           LQ258281         Negetive         R         R73400111         2016/06/21         16.25.32 PM         Indsayw         I	Facility Name: Address:	Immucor Bidg 700 Norcross	R&D 0 , GA	Assay ID: 3	_Cell	Instrument SN Software Versi Assay Version Report Version	: 503009087 ion: 1.6.11.1 : 4.00 05No n: 1.02	73 2.11 SP w12 3.00 30	Apr14
Visitive         PR73400111         2016/06/21         18:25:32 PM         Indisayw         Indisayw           Indicative         Indicative	Sample ID:		Interpretation:	Reference	Plate ID:	Date Read:	Time Read:	Operator ID:	Flag
y:       *INV* = invalid Interpretation No_int = No Integretation C/CIF # No Integretation       Figs:       2 = Warring 3 = Error	LQ258281		Negative		R73400111	2016/06/21	18:25:32 PM	Indsayw	
9. "INV" = Invalid Interpretation No_Int = No Interpretation CITIF all No Interpretation CITIF all No Interpretation CITIF all No Interpretation									
9.         "INV" = invalid Interpretation No_Int = No Integretation CITIF all Control Fillure         Fillings         2         = Warring 3         = Error									
9. "NV" = invalid Interpotation No_int = No Integretation CIT Falle Control Files									
y:     *INV* = invalid Interpretation No_int = No Integretation C/CIT Falle Control Files     Files:     2 = Warring 3 = Error									
y:         *INV* = Invalid Interpretation No_int = No Interpretation Cit F = No Interpretation Cit F = No Interpretation         Figs: 2 = Warring 3 = Error									
9. "NV" = Invalid Interpretation No_Int = No Integretation No_Int = No Integretation CIT Fall = Control Files									
9. "NV" = invalid Interpotation No_int = No Integretation Cit Fail = Control = Flags: 2 = Warring 3 = Error									
y: "NV" = invalid Interpretation No_Int = No Interpretation CIT F all No Interpretation									
y: "NV" = Invalid Interpretation No_int = No Interpretation Crit Fall = Control Files									
9. "INV" = Invalid Interpretation No_Int = No Interpretation CIT Fall = Control Files									
9. "NV" = invalid Interpretation No_int = No Integretation No_int = No Integretation Cit Fall = Control Files									
y: "NV" = invalid Interpretation No_int = No Interpretation No_int = No Interpretation CITIF = No Interpretation									
y: "NV" = Invalid Interpretation No_int = No Interpretation Crit Fall = Control Files Files: 2 = Warring 3 = Error									
y: "NV" = Invalid Interpretation No_Int = No Interpretation CIT Fall = Control Files									
y: "NV" = Invalid Interpretation No_Int = No Integretation Crit Field Control Fielder									
y: "NV" = Invalid Interpretation Flags: 2 = Warring No_Int = No Interpretation 3 = Error									
y: "INV" = Invalid Interpretation No_Int = No Interpretation Crit Fall = Control Fallure									
y: "INV" = Invalid Interpretation No_Int = No Interpretation Crit Fall = Control Fallure									
y: "INV" = Invalid Interpretation No_Int = No Interpretation Crit Field Control Field Field									
y: "INV" = Invalid Interpretation Flags: 2 = Warring No_Int = No Interpretation 3 = Error Crti Falve									
y: "INV" = Invalid Interpretation Flags: 2 = Warring No_Int = No Interpretation 3 = Error Chi Fail = Control Failure									
y: "INV" = Invalid Interpretation Flags: 2 = Warring No_Int = No Interpretation 3 = Error Cit Fall = Control Fallwe									
y: "INV" = Invalid Interpretation Flags: 2 = Warring - No_Int = No Interpretation 3 = Error Crt Falue									
No_Int = No Interpretation 3 = Error Cht Fail = Control Failure	Key: 1	NV" = In	velid Interpretatio	n			Flags	2 = Warring	
Ctri Fail = Control Failure	N	io_int = Na	Interpretation					3 = Emor	
		H Enil - C	ontrol Enilum						
	C	arran - G	onoor Panare						
E09/00 16-02-07 DM	C	arran - G	ond Panale	16-00	-77 DM				

## Single Assay Sample Based Report Columns

The table below describes the Sample Report columns.

Column	Description
Sample ID	The patient identification barcode information as it appears on the tube.
Interpretatio n	The most current interpretation for the sample. ABO and Rh columns report separate interpretations for the blood type.
Reference	If the current interpretation is from an edit, the report displays 'Manual Edit' in this column. If the current interpretation is from the result of a validate method, the report displays 'VALIDATE' in this column.
Plate ID	The plate ID from the location from which the sample results were generated.
Date Read	The date the NEO Iris reader read the plate.
Time Read	The time the NEO Iris reader read the plate.
Operator ID	The ID of the operator logged in when the sample results were generated. If the results were edited or run with a validate method, the report displays the ID of the user who was logged in at the time of the edit or running of the validate method.
Flag	This processing flag table cell is either left blank or is populated with a 2 or a 3. The key states that 2 indicates a <b>Warning</b> and 3 indicates an <b>Error</b> . A blank cell indicates no associated flags.

## Current Reports

#### **Current Report Types**

Current reports are a variation of sample based reports and are accessed as such. Current reports cannot be accessed through the plate based method of printing reports. Only one sample is listed per page along with the reagents used (not including the control information). Plate based reports are the only method to print control results. The result printed on the current report is the most current result. The three (3) broad categories of current reports are:

- Non-antibody identification reports
- Antibody identification reports
- Group and screen reports

#### Non-antibody Identification Current Report Details

This sample based report lists one sample per page, along with the reagents used (not including the plate control information). The assay name is displayed in the report header.

Facility Name: Im Address: Bi No	imucor R&D dg 7000 proross, GA					Assay Opera Date F Time I	r ID: itor ID: Read: Read:	ABORH lindsayw 2016/06/21 17:37:58 PN	л		Instrument SN Software Versi Assay Version Report Version	: 50 ion: 1.6 : 3.0 : 1.0	30090873 i.11.1 i2i 28Jan16 i9	2.11 SP1 2.03i 04Feb16
Plate ID: U/	11139582	-							44.0-5	0.0-1	D-f	Fine		
Sample ID:	ABO	Rh	nitheo	Anti-A	Anti-B	Anti-D4	Anti-D5	Mono Cen	AI-Cell	A B-Cell	Reference	Flag		
A1-Cell (10ml) Anti-A (10ml) B-Cell (10ml)	11 10	1083 1746 3083	2016/07/08 2017/09/01 2016/07/08											
Anti-A (10ml)	10	1746	2017/09/01	_										
Anti-B (10ml)	20	3620	2010/07/08	-										
Anti-D series 4 (1	10ml) 50	4715	2018/03/03	-										
Anti-D series 5 (1	10ml) 50	5741	2017/08/11											
Key: "INV" NTD -	= Invalid Ir = No Type = Negative	nterpretatio Determine Result	n 1-4 = G d ? = E X = In	raded P quivocal valid Re	ositive Re Result sult	esult						Flags:	2 = W 3 = Er	faming tror
2016/06/22			16:03:3	9 PM										Page: 1

## **Antibody Identification Report Details**

The antibody identification reports have a unique format to better align with Masterlists. The plate ID that sample was run with is located on the left side of the report directly above the results grid. Only one sample is listed per page along with the reagents used.

IMMUC	OR.	Ab_ID C	Current Res	ult Report	Z	EO t IRIS	
Facility Name: Immuc Address: Bidg 7 Norcro	or R&D 000 ss, GA	Assay ID: Operator ID: Date Read: Time Read:	Ab_ID lindsayw 2016/06/21 17:48:30 PM	Instrument SN: Software Version Assay Version: Report Version:	5030090873 : 1.6.11.1 4.00 05Nov12 1.14	2.11 SP1 1.05 30Apr14	
Lot Number: ID301 Plate ID: ID3010	)2556						
Reagent Name: Cap-R Pos Ctrl (10m	Lot Number: () 244330	Expiry Date: 2017/01/15	Sample	Number	LQ258281		
LISS (10ml)	211750	2017/03/29	Interpret	ation	Negative		
Cap-R Ind Cells (10n	N) 221650	2016/07/07	Cell 1		-		
Cap-R Neg Ctrl (10m	1) 245190	2017/01/14	Cell 2		_		
			Cell 3		-		
Attention:			Cell 4		-		
This report shows same	ple reactivity base	ed on a specific	Coll 5		-		
lot of test plates. It is i	mportant that the	Master List for	Cell 6		-		
that plate lot is used w	nen interpreting th	rese results. Use	Cell 6		-		
misidentification of the	antibody present	in the second se	Cell /		-		
			Cell 8		-		
The plate lot used for t	his test is:						
Lot	No: ID30	1	Cell 9				
		-	Cell 10		-		
			Cell 11		-		
			Cell 12		-		
			Cell 13		-		
			Cell 14		-		
			Positive	Control	4		
			Negative	e Control	-		
			Referen	се			
			Flag				
Key: "INV* =   No_Int =	Invalid Interpretation No Interpretation Negative Result	on ? = X = Ctrl Fail :	Equivocal Result Invalid Result = PC < 3+ or NC >	Flags:	2 = Warning 3 = Error		
- 1-4 =	Graded Positive R	lesult					

## Antibody Identification Report Columns and Features

The table below describes the antibody identification report columns and features.

Column or feature	Description
Reagent Name:, Lot Number: and Expiry Date	A reagent details grid lists the name lot number and expiry date used in the assay.
Sample Number	The patient identification barcode information as it appears on the tube.
Interpretation	The most current interpretation for this sample, run with this assay.
Result text	The result text grid lists each test well and the result that was assigned by the software. Each test well is described, e.g. Cell 13.
Reference	If the current interpretation is from an edit, the report displays 'Manual Edit'. If the current interpretation from the result of a validate method, the report displays 'VALIDATE'.
Flag	This processing flag table cell is either left blank or is populated with a 2 or a 3. The key states that 2 indicates a <b>Warning</b> and 3 indicates an <b>Error</b> . A blank cell indicates no associated flags.
Ctrl Fail = PC < 3+ or NC > 0	Some current reports have criteria set for qualifying the control results. These criteria are listed in the footer key. If actual results fall outside of the criteria, then the controls are deemed to have failed and are invalid.
Attention:	As an example, red blood cell antibody identification reports include an attention that gives guidance on the use of the report. It states: "This report shows sample reactivity based on a specific lot of test plates. It is important that the Master List for that plate lot is used
	when interpreting these results. Use of the Master List of a different lot could result in misidentification of the antibody present."
Plate lot number	The plate lot number is reiterated in large bold text to stress this information.

#### Group and Screen Sample Based Report Details

The Group and Screen report combines sample results from two assays onto one page. The sample ID for which the instrument generates results is located on the left side of the group and screen report, directly above the results grid. This report lists results for sample runs with the ABORH or the ReflexABO assay combined with a screen assay (Pool\_Cell, 2\_Cell, or 3\_Cell). If results exist for only one of the assays, the other grid table cells remain blank (refer to the example report shown below). This report lists one sample per page, along with the reagents used (not including the control information) for both assays.

										,							IRIS		
Facility Name:	Immucor	R&D										Instrum	ent SN:		5030090	873			
Address:	Bidg 7000	1										Softwar	re Version:		1.6.11.1		2.11 SP	1	
	Norcross,	GA										Assay	Version (ABC	PRH):	3.021 283	anito	2.03104	April	
												Assay	Version (5_C	en):	4.00 054	40412	3.00 30	Mpi 14	
Sample ID:	LQ258281											кероп	version:		1.02				
Assay Name:	ABO	Rh		Anti-A	Anti-B	Anti-D4	Anti-D5	Mono Ctrl	A1-Cell	B-Cell	Refere	ince	Plate ID:	0	ate Read:	Ti	me Read:	Operator ID:	Fla
ABORH	0	Positive				4	4		4	4			UA111395	82 2	016/06/21	17	:37:58 PM	lindsayw	
Assay Name:	Reage	ent Name:		Lot Nu	mber:	Lot Expin	y:												
Assay Name:	Reage	ent Name:		Lot Nu	mber:	Lot Expin	y:												
Assay Name: ABORH	Mono	ctrl (10ml)		49224	mber:	Lot Expiry 2018/03/0	y: D1												
Assay Name: ABORH ABORH	Mono A1-Ce	ent Name: ctrl (10ml) ell (10ml)		49224 111083	mber:	Lot Expiry 2018/03/0 2016/07/0	y: D1 D8												
Assay Name: ABORH ABORH ABORH	A1-Ce Anti-A	ent Name: ctrl (10ml) ell (10ml) (10ml)		49224 111083 101746	mber:	Lot Expiry 2018/03/0 2016/07/0 2017/09/0	y: 01 08 01												
Assay Name: ABORH ABORH ABORH ABORH	Anti-A B-Cell	ent Name: ctrl (10ml) ill (10ml) (10ml) (10ml) (10ml)		Lot Nu 49224 111083 101746 113083	mber:	Lot Expiry 2018/03/0 2016/07/0 2017/09/0 2016/07/0 2017/10/2	y: 01 08 01 08												
Assay Name: ABORH ABORH ABORH ABORH ABORH	Anti-A Anti-A Anti-A	ent Name: ctrl (10ml) ill (10ml) (10ml) (10ml) (10ml) series 4 (1)	Omili	Lot Nu 49224 111083 101746 113083 203620 504716	mber:	Lot Expiry 2018/03/0 2016/07/0 2017/09/0 2016/07/0 2016/07/0 2018/03/0	y: 01 08 01 08 20												
Assay Name: ABORH ABORH ABORH ABORH ABORH ABORH	Anti-A Anti-A B-Cell Anti-B Anti-D Anti-D	ent Name: ctrl (10ml) (10ml) (10ml) (10ml) series 4 (1) series 5 (1)	Oml)	Lot Nu 49224 111083 101746 113083 203620 504710 50574	mber:   	Lot Expiry 2018/03/0 2016/07/0 2017/09/0 2016/07/0 2017/10/2 2018/03/0 2017/08/1	y: 01 08 01 03 20 03												
Assay Name: ABORH ABORH ABORH ABORH ABORH ABORH 3 Cell	Reage Mono A1-Ce Anti-A B-Cell Anti-B Anti-D Anti-D LISS	ent Name: ctrl (10ml) (10ml) (10ml) (10ml) series 4 (1) series 5 (1) 10ml)	Oml) Oml)	Lot Nu 49224 111083 101746 113083 203620 504713 50574 211750	mber:   	Lot Expiry 2018/03/0 2016/07/0 2017/09/0 2016/07/0 2016/03/0 2018/03/0 2017/08/1 2017/03/2	y: 01 08 01 08 20 03 11 29												
Assay Name: ABORH ABORH ABORH ABORH ABORH ABORH 3_Cell 3 Cell	Reage Mono A1-Ce Anti-A B-Cell Anti-D Anti-D LISS ( Cap-R	ent Name: ctrl (10ml) (10ml) (10ml) (10ml) series 4 (1) series 5 (1) 10ml) Lind Cells (	0ml) 0ml) 10ml)	Lot Nu 49224 111083 101746 113083 203620 504711 50574 211750 221650	I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I	Lot Expiry 2018/03/0 2016/07/0 2017/09/0 2017/09/0 2017/10/2 2018/03/0 2017/08/1 2017/08/1 2017/08/2 2016/07/0	y: 01 08 01 08 20 03 11 29 07												
Assay Name: ABORH ABORH ABORH ABORH ABORH ABORH 3_Cell 3_Cell	Reage Mono A1-Ce Anti-A B-Cell Anti-B Anti-D Anti-D LISS ( Cap-R	ent Name: ctrl (10ml) (10ml) (10ml) (10ml) series 4 (1) series 5 (1) 10ml) t Ind Cells (	Oml) Oml) 10ml)	Lot Nu 49224 111083 101746 113083 203620 504710 50574 211750 221650	I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I           I         I	Lot Expiry 2018/03/0 2016/07/0 2016/07/0 2016/07/0 2017/10/2 2017/08/1 2017/03/2 2016/07/0 2016/07/0	y: 001 008 001 008 000 003 003 011 11 229 007												
Assay Name: ABORH ABORH ABORH ABORH ABORH ABORH 3_Cell 3_Cell	Reage Mono A1-Ce Anti-A B-Cell Anti-D Anti-D LISS ( Cap-F	ent Name: ctrl (10ml) (10ml) (10ml) (10ml) (10ml) series 4 (1 series 5 (1 (10ml) t Ind Cells (	0ml) 0ml) 10ml)	Lot Nu 49224 111083 101746 113083 203620 504718 50574 211750 221650	I         I           I         I	Lot Expin 2018/03/0 2016/07/0 2016/07/0 2016/07/0 2017/10/2 2018/03/0 2017/08/1 2017/03/2 2016/07/0	y: 001 008 001 008 000 003 003 003 007 007												
Assay Name: ABORH ABORH ABORH ABORH ABORH ABORH 3_Cell 3_Cell 3_Cell	Reage       Mono       Anti-A       B-Cell       Anti-D       Anti-D       Anti-D       LISS I       Cap-Fr       V*	ent Name: ctrl (10ml) ill (10ml) (10ml) (10ml) (10ml) i series 4 (1 i series 5 (1 (10ml) t Ind Cells ( lid Interpret	0ml) 0ml) 10ml) ation	Lot Nu 49224 111083 101744 113083 203620 504712 50574 211750 221650 <b>1-4</b> =	mber:	Lot Expin 2018/03/0 2016/07/0 2016/07/0 2016/07/0 2017/10/2 2018/03/0 2017/08/1 2017/03/2 2016/07/0 Positive I	y: 01 08 01 08 01 08 020 03 03 11 29 07 Result					Flags:	2 = Wan	ning					
Assay Name: ABORH ABORH ABORH ABORH ABORH ABORH 3_Cell 3_Cell 3_Cell 	V* = Inva Int = No I	ent Name: ctrl (10ml) ill (10ml) (10ml) (10ml) (10ml) series 4 (1 series 5 (1 10ml) t Ind Cells ( lid Interpret nterpretatio	0ml) 0ml) 10ml) ation	Lot Nu 49224 111083 101744 113083 203620 504713 50574 221650 221650 <b>1-4</b> = <b>?</b> =	I         I           I         I	Lot Expin 2018/03/0 2016/07/0 2016/07/0 2016/07/0 2017/09/2 2016/07/0 2017/03/2 2016/07/0 Positive I al Result	y: 01 08 01 08 00 03 03 11 29 07 Result					Flags:	2 = Wan 3 = Error	ning					
Assay Name: ABORH ABORH ABORH ABORH ABORH ABORH ABORH ABORH 3_Cell 3_Cell 3_Cell No. NT	Reage       Mono       A1-Ce       Anti-A       B-Cell       Anti-D       LiSS (       Cap-F       V*       Int = No I       D       No I	ent Name: ctrl (10ml) ill (10ml) (10ml) (10ml) (10ml) iseries 4 (1 series 5 (1 10ml) lind Cells ( lid Interpret nterpretatio fype Determ	Oml) Oml) 10ml) ation n	Lot Nu 49224 111083 101744 113083 203620 504710 50574 211750 221650 <b>1-4=</b> <b>?</b> = <b>X</b> =	Graded Fundamental Fundamentar	Lot Expin 2018/03/0 2016/07/0 2016/07/0 2016/07/0 2016/07/0 2017/03/2 2016/07/0 2017/03/2 2016/07/0 Positive I al Result Result	y: 11 12 13 14 15 16 17 17 17 17 17 17 17 17 17 17					Flags:	2 = Wan 3 = Erro	ning					

#### Group and Screen Sample Based Report Columns

The table below describes the Group and Screen Report columns.

Column	Description
Assay Name	The assay name for the results generated.
Interpretatio n	The most current interpretation for this sample run with this assay. ABO and Rh columns report separate interpretations for the blood type.
Result Text	The result text column(s) represent each test well (with a description of the well, such as Anti-A) and the result that the system used to determine the interpretation.

Column	Description
Reference	If the current interpretation is from an edit, the report displays 'Manual Edit' in this column. If the current interpretation is from the result of a validate method, the report displays 'VALIDATE' in this column.
Plate ID	The plate ID from the location from which the sample results were generated.
Date Read	The date the NEO Iris reader read the plate.
Time Read	The time the NEO Iris reader read the plate.
Operator ID	The ID of the operator logged in when the sample results were generated. If the results were edited or run with a validate method, the report displays the ID of the user who was logged in at the time of the edit or running of the validate method.
Flag	This processing flag table cell is either left blank or is populated with a 2 or a 3. The key states that 2 indicates a <b>Warning</b> and 3 indicates an <b>Error</b> . A blank cell indicates no associated flags.

## Quality Control Reports

#### **Report Details**

Quality Control (QC) reports are only available as plate based reports and not through sample based report printing methods. The plate ID for which the quality control (QC) results are generated is located on the left side of the QC report, directly above the result grid. The result grid identifies the reagents in each test well, the result generated with this run, and the expected result.

The reagent details grid is located directly below the result grid. This report lists the reagents used for this QC with their lot number, expiry date and reagent barcode information. The **QC Result** column lists either 'Qualified' or 'Failed' as a summary, depending on the outcome of the individual results and if those individual results match the expected results.

There are multiple QC assays dependent upon the selected reagents being used for different blood group purposes. Refer to **Attachment 1 for NEO Iris Operator Manual** for a detailed list of the QC assays and their applications.

In cases where multiple vials of the same lot of reagent are loaded and used in single quality control test assays such as QCTEST, then the report supports the documentation of multiple vials.

The entire barcode information of the used reagent vials is documented on all quality control reports along with the lot number and expiry date of those same vials.

Below are two (2) examples of QC reports.

Facility Name: Immucor R&D Address: Bldg 7000 Norcross, GA		Assay II Operato Date Re Time Re	): QCTES r ID: lindsay ad: 2016/0 pad: 17:09:4	ST I W S 6/21 A 42 PM F	nstrument SN: Software Version: Assay Version: Report Version:	5030090873 1.6.11.1 3.04i 11Apr16 1.11	2.11 SP1 1.02 30Apr1
Plate ID: UA11139582					OC Resul	÷-	T
Lot Number: UA111 Expiry Date: 2030/12/30					Qualified		
Well	Result	Expec	ted Result				
A1: Anti-A & A1-Cells	+	+					
B1: Anti-A & B-Cells							
C1: Anti-B & B-Cells	+	+					
D1: Anti-B & A1-Cells		-					
E1: Anti-D4 & corQC Extend STD	+	+					
F1: Anti-D4 & A1-Cells	-						
G1: Anti-D5 & corQC Extend STD	+	+					
H1: Anti-D5 & B-Cells		-					
A1-Cell (10ml)	1110	83	2016/07/08	1111901608303	323687		
Mono ctri (10mi)	4922	:41	2018/03/01	492060182410	411186		
B-Cell (10ml)	1130	.83	2016/07/08	113190160830	315800		
Anti-B (10ml)	2036	120	2017/10/20	203293176201	672615		
Anti-D series 4 (10ml)	5047	15	2018/03/03	504062187151	529545		
Anti-D series 5 (10ml)	5057	741	2010/07/08	605222177410	190794		
Prior January (Torin)	3037	41	2017/00/11	505225177410	370396		

Facility Name: Immus Address: Bidg 7 Nororo Plate ID: R7340 Lot Number: R734 Expiry Date: 2016/0 Well A1: Expected Positiv B1: Expected Positiv C1: Expected Positiv E1: Expected Positiv E1: Expected Positiv E1: Expected Negati F1: Expected Negati G1: Expected Negati	cor R&D r000 coss, GA 00111 07/15 Result re + re +	Assay ID: Operator I Date Rea: Time Rea Time Rea Expected Result + + + + - - -	QC3_Ce D: lindsayw ± 2016/06/4 d: 17:52:52	QC R	Instrument SN: Software Version: Assay Version: Report Version: Result: ified	5030090873 1.6.11.1 2.11 SP1 4.00 23Nov12 3.01 18Dec1 3.00
Address: Bidg 7 Nororo Plate ID: R7340 Lot Number: R734 Expiry Date: 2016/0 Well A1: Expected Positiv B1: Expected Positiv C1: Expected Positiv E1: Expected Positiv E1: Expected Negati F1: Expected Negati F1: Expected Negati H1: Expected Positiv	Result           re         +           re         +           re         +           ive         -	Operator I Date Rea: Time Rea: Time Rea: Expected Result + + + - - - +	D: lindsayw d: 2016/06/ d: 17:52:52	QC R Quali	Software Version: Assay Version: Report Version: Result: ified	1.6.11.1 2.11 SP1 4.00 23Nov12 3.01 18Dec1 3.00
Noraco Plate ID: R7340 Lot Number: R734 Expiry Date: 2016/0 Well A1: Expected Positiv B1: Expected Positiv D1: Expected Positiv E1: Expected Positiv E1: Expected Negati F1: Expected Negati H1: Expected Positiv	Result           77/15           Result           re           +           re           re      re	Date Real Time Real Expected Result + + + + - - - - +	d: 2016/06/ d: 17:52:50	QC R Quali	Assay Version: Report Version: Result: ified	4.00 23Nov12 3.01 18Dect 3.00
Plate ID: R7340 Lot Number: R734 Expiry Date: 2016/C Well A1: Expected Positiv B1: Expected Positiv C1: Expected Positiv C1: Expected Positiv E1: Expected Negati F1: Expected Negati F1: Expected Negati H1: Expected Positiv	00111 07/15 Result re + re	Time Real Expected Result + + + - - - +	d: 17:52:5	QC R Quali	Report Version: Result: ified	3.00
Plate ID: R7340 Lot Number: R734 Expiry Date: 2015/0 Well A1: Expected Positiv B1: Expected Positiv D1: Expected Positiv D1: Expected Positiv E1: Expected Negati F1: Expected Negati F1: Expected Negati H1: Expected Positiv	Result           re         +           re         +           re         +           ve         +           vive         -           live         -           vive         -           vive         -           vive         -           vive         -           vive         -	Expected Result +		QC R Quali	Result: ified	
Lot Number: R734 Expiry Date: 2016/C Well A1: Expected Positiv C1: Expected Positiv C1: Expected Positiv D1: Expected Positiv E1: Expected Negati F1: Expected Negati H1: Expected Positiv	Result           re         +           re         +           re         +           ve         +           vive         -           vive         -           veve         +	Expected Result + + + - - - +		Quali	ified	
Vell A1: Expected Positiv B1: Expected Positiv C1: Expected Positiv D1: Expected Positiv D1: Expected Positiv E1: Expected Negati G1: Expected Negati H1: Expected Positiv	Result           re         +           re         +           ve         +           vive         -           vive         -           vive         -           vee         +	Expected Result + + + + - - - - + +		Quali	ified	
Well A1: Expected Positiv B1: Expected Positiv C1: Expected Positiv D1: Expected Positiv E1: Expected Negati F1: Expected Negati H1: Expected Positiv	Result           re         +           re         +           ve         +           vive         +           vive         -           vive         -           vive         -           vive         -           vive         -           vive         -	Expected Result + + + + + + + + + + + + + +				
A1: Expected Positiv B1: Expected Positiv C1: Expected Positiv D1: Expected Positiv E1: Expected Negati G1: Expected Negati H1: Expected Negati H1: Expected Positiv	ve     +       ve     +       ve     +       ve     -       ive     -       ive     -       ve     +	* * * * * * * * * * * * * * * * * * *				
B1: Expected Positiv C1: Expected Positiv D1: Expected Positiv E1: Expected Negati G1: Expected Negati H1: Expected Positiv	ve + ve + ive - ive - ive - tive - ve +	+ + - - - +				
C1: Expected Positiv D1: Expected Positiv E1: Expected Negati F1: Expected Negati G1: Expected Negati H1: Expected Positiv	ve + ve + ive - ive - ive - ve +	* * - - *				
D1: Expected Positiv E1: Expected Negati F1: Expected Negati G1: Expected Negati H1: Expected Positiv	ve + ive - ive - tive - ve +	+ - - +				
E1: Expected Negati F1: Expected Negati G1: Expected Negati H1: Expected Positiv	ive - ive - tive - ve +	• • •				
F1: Expected Negati G1: Expected Negati H1: Expected Positiv	ive - tive - ve +	- -				
G1: Expected Negati H1: Expected Positiv	tive - ve +	+				
H1: Expected Positiv	ve +	+				
Reagent Name		Lot Number	Expiry Date	Barcode		
Cap-R Ind Cells (10r	mi)	221650	2016/07/07	2211891665	500565619	_
LISS (10ml)		211750	2017/03/29	2110881775	501027406	
Cap-R Pos Ctrl (10m	nl)	244330	2017/01/15	2440151733	300286286	
Cap-R Neg Ctrl (10n	ni)	245190	2017/01/14	2450141719	900241774	-

## Reagent Reports

There are two (2) versions of reports that list reagent details. One version of reagent report is derived from plate based reports and is not accessible through the sample based methods of reporting. This version uses the *Reporting* dialog. The other version of reagent report is derived from both sample and plate based methods of reporting using the **Details** button.

## Using the Details Button

You can access one version of **Reagents** reports by selecting the *Reagents* tab of the *Test Details* dialog, via the **Details** button, for a chosen sample or plate item on the main *Results* screen.

	Ē	Actions -	1
		0	l
Details		Details	l

well	Reagent	Lot Number	Qualified	Expiration	Last QC	QC Expires	Barcode	Edited	Edited B
A02	Anti-A	101711	1	2011/03/23	2010/02/25 08:01:47	2010/02/26 08:01:32	101082117110949625		
302	Anti-B	203281	1	2011/05/07	2010/02/25 08:01:47	2010/02/26 08:01:32	203127112819824198		
02	Anti-D series 4	504732	1	2011/03/07	2010/02/25 08:01:47	2010/02/26 08:01:32	504066117321675430		
02	Anti-D series 5	505571	1	2011/05/07	2010/02/25 08:01:47	2010/02/26 08:01:32	505127115711037806		
02	Mono ctrl	492051	1	2011/08/06	2010/02/25 08:01:47	2010/02/26 08:01:32	492218110511853900		
02	A1-Cell	111737	1	2010/02/26	2010/02/25 08:01:47	2010/02/26 08:01:32	111057107370326511		
302	B-Cell	113737	1	2010/02/26	2010/02/25 08:01:47	2010/02/26 08:01:32	113057107370185817		

Press the Print button to print the Reagents report.



An example of a **Reagents** report derived from the *Reagents* tab is shown below. This version of reagent report lists the reagents used with a given sample or plate. The sample number field displays the sample identification that is currently selected or is defaulted to the top sample in any sample list if a specific sample is not selected, such as in the case of plate based reports. In the case of quality control assays, the sample number field displays the barcoded identification information of one of the reagents undergoing quality control on that quality control plate.

IRIS		<b>Re</b> Database IMM BL	Reagents Database: Galileo_SQL IMMUCOR R&D BLDG 7000					
te Number: UA111399 say Name: FWDABO mple Number: W05331	90 _AB 6502593							
te Number: UA111399 say Name: FWDABO mple Number: W05331 Reagent Name	90 _AB 6502593 Lot No.	Expiration	Barcode	Well(s)				
te Number: UA111399 say Name: FWDABO mple Number: W05331 Reagent Name Anti-A (10m)	90 _AB 6502593 Lot No. 101725 200590	Expiration 2017-03-17 2017-03-24	Barcode 101076177251060795 203083175001976130	Well(s) A01A06 B01 B06				
tte Number: UA111399 say Name: FWDABO mple Number: W05331 Reagent Name Anti-A (10ml) Anti-B (10ml) Anti-B (10ml)	90 _AB 6502593 Lot No. 101725 203590 301430	Expiration 2017-03-17 2017-03-24 2017-04-15	Barcode 101076177251060795 203083175901976130 301105174300885797	Well(s) A01A06 B01B06 C01C06				
te Number: UA111399 say Name: FWDABO, mple Number: W05331 Reagent Name Anti-A (10m) Anti-B (10m) Anti-B (10m) Anti-B (10m)	90 AB 6502593 Lot No. 101725 203590 301430 504971	Expiration 2017-03-17 2017-03-24 2017-04-15 2017-04-07	Barcode 101076177251060795 203083175901976130 301105174300885797 504097179711744147	Well(s) A01A06 B01B06 C01C06 D01D06				

Refer to the **Printing Reagents Report from Test Details** section of this chapter for details of this method of printing.

## Using the Reporting Dialog

The second version of a report listing reagents (the **Reagent Report**) is available using the *Reporting* dialog, as described in the **Printing Reports** section of this chapter. This version of **Reagent Report** lists the reagents used on a given plate on a sample-by-sample basis for all of the sample identifications included on that plate. This version of the report is not accessible through the sample based methods of reporting.

MMUC	OR.			Reagen	t Report			NEO + IRIS	
Facility Name: Imm Address: Bidg Norr Plate ID: ID30 Lot Number: ID30	ucor R&D 17000 27055, GA 1102556 11			Assay ID: Date Read: Time Read:	Ab_ID 2016/06/21 17:48:30 PM		Instrument SN: 50300 Software Version: 1.6.11 Report Version: 1.07	990873 1.1 2	.11 SP1
Sample ID:	Respect Name:	Lot Number	Evniry Date:	Respect Name	Lot Numbe	Expiry Date:	Respert Name:	Lot Number	Expiry Date:
LQ258281	Cap-R Ind Cells (10ml)	221650	2016/07/07	LISS (10ml)	211750	2017/03/29	Cap-R Pos Ctrl (10ml)	244330	2017/01/15
	Cap-R Neg Ctrl (10ml)	245190	2017/01/14						

Refer to the Printing Reports section of this chapter for details of this method of printing.

## Accessing Plate Based Reports

You can access plate based reports by setting the results view to plate, selecting the plate entries to be reported, and pressing the **Print** button in the upper right corner.



View	Data Selec Plate : Sample :	All All	Assay : All From : 2009/09	Appr     Appr     Appr     All     C No     C No     C No     C No     C No	oval Status Tests of Reviewed oproved Onl of App <u>r</u> oved		t Types Types iusual grmal	Actions Details	Approve Print	Export
« Date	Conductri	Plate Read Date	Plate ID	Assav Name	Tests F	ans Edit	d Approved	Edit		Refresh
2009/09/17	7 10:06:55	2009/09/17 10:06:	35 LIA5364570	mAgScr37	91	ago   care	a rippieres	Landoucord	Plate has interpreted	without
2009/09/17	7 09:19:19	2009/09/17 09:18:4	48 SC12102915	5 WEAK D	46				Plate has interpreted	without
2009/09/17	7 08:10:39	2009/09/17 08:10:2	24 UA5364569	OCTEST	8				Plate has interpreted	without
2009/09/16	5 14:27:16	2009/09/15 09:50:	25 UA5415035	ABORH	2	,			Plate has interpreted	without
2009/09/19	514:26:03	2009/09/15 14:25:	53 DP03501007	7 ExtendDP	1				Plate failed criteria.	lease ch
2009/09/15	5 13:00:45	2009/09/15 13:00:	34 SC12102671	WEAK D	1				Plate has interpreted	without
2009/09/19	5 12:56:53	2009/09/15 12:56:	44 UA5364519	ABORH	1				Plate has interpreted	without
2009/09/15	5 12:33:38	2009/09/15 12:33:1	15 UA5364519	OCTEST	8				Plate has interpreted	without
2009/09/19	5 09:50:38	2009/09/15 09:50:	25 UA5415035	ABORH	2				Plate has interpreted	without
2009/09/14	14:00:43	2009/09/14 14:00:2	26 UA5415034	FWD ABORH	4	1			Plate has interpreted	without
2009/09/14	13:54:45	2009/09/14 13:54:	35 UA5415085	FWD ABORH	12				Plate has interpreted	without
2009/09/14	13:50:03	2009/09/14 13:49:4	47 UA5415086	FWD ABORH	12				Plate has interpreted	without
2009/09/14	13:45:53	2009/09/14 13:45:4	40 UA5415083	FWD ABORH	12				Plate has interpreted	without
2009/09/14	13:40:42	2009/09/14 13:40:2	21 UA5415084	FWD_ABORH	12				Plate has interpreted	without
2009/09/14	13:34:09	2009/09/14 13:33:5	54 UA5415090	FWD ABORH	12				Plate has interpreted	without
•					1.11777.1					×
Database	Name: G	alileo-PC/Galileo	_SQL							OK



Note: You must select at least one plate entry before pressing the Print button.

You can also access plate based reports by selecting the **Print** button in the *Plate Views* tab of the *Test Details* dialog of a plate. You can access *Test Details* by pressing the **Details** button.



ample No	Edited ABO	Rb	Well	View	
R132522	В	Positive	A02H02		1
R132493	0	Positive	A03H03	T D Reaction	1
R132559	0	Positive	A04H04	+ / - Grade Strength	I
R132521	А	Positive	A05H05		
R132488	0	Positive	A06H06	1 2 3 4 5 6 7 8 9	10
R132554	0	Positive	A07H07		$\square$
R132547	А	Positive	A08H08		X
R132532	А	Positive	A09H09		+)
R132517	В	Positive	A10H10		$\overline{\mathbf{n}}$
R132494	А	Positive	A11H11		<b>T</b>
R132527	0	Negative	A12H12	⋼∩⋳⋳⋳⋳⋳	<del>A</del>
					*
					-)
					Š.
					J
					$\neg$
					$\asymp$
					)(
				1000000000	-
				Approve Export Print	

## Accessing Sample Based Reports

You can access sample based reports by setting the results view to sample, selecting the sample entries to be reported, and pressing the **Print** button in the upper right corner.

and the second	-	
	Print	

View	Data Selec Plate : Sample : Donation :	All · · · · · · · · · · · · · · · · · ·	Assay From To	: All : 2009/09/08 : 2009/09/17	•	Approval Status All Tests Not Reviewer Approved Or Not Approver Current Only	Result Typ All Typ Unusua Normal Other Eilte	pes es al ers	Actions Details App Edit	prove Print <u>Y</u> oid	Export 2 Refresh
Date »		Plate Read Date	S	ample ID		Assay Name	Result 1	Result 2	Donation ID	Flags Edited	Approvec 🔺
2009/09/10	11:55:11	2009/09/10 11:54:	36 R	133852		Pool_Cell	Negative	10		82 - 22.	
2009/09/10	11:55:11	2009/09/10 11:54:	36 R	133812		Pool_Cell	Negative				
2009/09/10	11:55:11	2009/09/10 11:54:	36 R	133807		Pool_Cell	Negative				
2009/09/10	11:55:11	2009/09/10 11:54:	36 R	133814		Pool_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133824		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133881		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133839		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133888		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133828		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133904		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133893		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133841		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133906		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133905		3_Cell	Negative				
2009/09/10	12:00:39	2009/09/10 12:00:	15 R	133838		3_Cell	Negative				-
4											•
Database i	Name: G	alileo-PC/Galileo	)_SQ	L							ОК



Note: You must select at least one sample entry before pressing the Print button.

You can also access sample based reports, such as Group and Screen reports, by using the **Print** button in the *Test Overview* tab or the *Results* tab of the *Test Details* dialog.



**Note**: Group and Screen reports can only be accessed using the sample based report methods and not by using the plate based report methods.

You can access the *Test Details* dialog by pressing the **Details** button in the *Actions* area of the main *Results* screen. You must select one sample entry before pressing the **Details** button.



#### Chapter 8: NEO Iris Reports

Sample ID : LQ136083	Gender :	Instrument ID : 5030090012 Software Version : 1.5.8-1.40 SP5
Patient ID :	Date of Birth : 2010/03/01	DMS Assay : 1.04 06July05
est Details Assay :  ABORH est Date :  2010/02/25 18:33:33	Microplate No. : UA5682662 Well Locations : [402H02	Test Status Approved Reported Exported
Current Results ABO O Rh Positive	Original Results	Edit
Comment :		

Test Details dialog>Test Overview tab

ample No.	Edited	ABO	Rh	Well		Well	Reagent	Well	Reaction	Original	Rev
2136083		0	Positive	A02H02	[	A02	ANTI-A		9	$\overline{\bigcirc}$	-
						B02	ANTI-B		14	Õ	-
						C02	ANTI-D4	۲	92	4	4
						D02	ANTI-D5		91	4	4
						E02	RHCTRL	۲	12	$\odot$	-
						F02	A1-CELL	۲	69	2	2
						G02	B-CELL	۲	67	2	2
						* H02	Empty	۲	$\overline{m}$	0	Γ
						AB	0	¥	] [Te	st Status Approv	ed
						R	h Positive	*	]	Report	ed ed
						Comment	;] 	1		L .	su

Test Details dialog>Results tab

## Printing Reports

## Printing Reports Using the Print Button

You can print reports by doing the following:

Step	Action					
1	Sample Based Reporting Options:					
	If you would like to print one of the following reports:					
	Sample based report					
	Antibody Identification report					
	Group and Screen report					
	The Print button should be selected from one of the following screens:					
	Sample View Main Results screen					
	• Test Details $\rightarrow$ Test Overview tab					
	• Test Details $\rightarrow$ Results tab					
	Plate Based Reporting Options:					
	If you would like to print one of the following reports:					
	Plate based report					
	Reagent report					
	Quality Control report					
	The Print button should be selected from one of the following screens:					
	Plate View Main Results screen					
	• Test Details $\rightarrow$ Plate View tab					

Step	Action
2	The system displays the <i>Reporting</i> dialog with the check box for the selected entry already checked.
	Report Name : Search Type:
	Plate         Assay         User         Date           INU15204562         Ag_C RH2         Donna         2010/01/27 10:48:29           Interview         Interview         Interview         Interview           Interview         Interview         Interview         Interview
	From the Report Name drop-down list, select the appropriate report.
	Reporting     x       Report Name :    Cell Current Result Report       2_Cell Sample Report
	Sample No.       4. Cell Sample Report         AD. D. Gurrent Result Report       ABO. AB.2 Current Result Report         ABO.AB.2 Current Result Report       ABORH.2 Current Result Report         ABORH.2 Sample Report       ABORH.2 Current Result Report         ABORH.2 Sample Report       ABORH.2 Sample Report         ADC. AB 2. Sample Report       ABORH.2 Sample Report         ADC. CE Current Result Report       ABORH.2 Sample Report         AD.C. E. Sample Report       AD.C. CE Current Result Report         AD.C. E. Sample Report       AD.C. E. Sample Report         AD.C. E. Current Result Report       AD.C. E. Sample Report         AD.C. E. Sample Report       AD.C. E. Sample Report         AD. C. E. Current Result Report       AD.C. E. Sample Report         AD. C. E. Current Result Report       C. C. Current Result Report         AD. C. E. Current Result Report       C. C. C. Current Result Report         AD. C. E. Sample Report       C. C. C. Current Result Report         AD. J. Ell Current Result Report       C. C. C. Current Result Report         AD. J. Ell S
	You can only print one report type at a time; therefore it is important to select a <i>Report Name</i> that corresponds to the assay name of the selected plate(s) or sample(s).

	Action	
The sample results to          Sample No.       Donation N         ☑       LQ136086         The Reporting dialog is         Reporting dialog is         Reporting dialog is         Sort/Elminate Sample None         Elminate Result : None         Sample No.       Donation No.         Assay         ☑       LQ136086         ☑       I         ☑       I         ☑       I         ☑       I         ☑       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I         I       I	be printed will be indicated by the box checked in front	of the ID.
То	press	Button
Select all sample IDs	Select All	V the to tail V the to tail V the to the V the to the V the to the Select All
Cancel the selection	<b>Cancel</b> . This action closes the <i>Reporting</i> dialog.	
Print the report	<b>Report</b> from the <i>Reporting</i> dialog after selecting	3

## **Automatic Printing of Failed Plate Reports**

A printed report is automatically generated when an assay with plate controls yields a control failed criteria. The automatically printed report indicates that the plate failed, what Control caused the failed status, and the list of the samples that were on the plate.



**Note**: The value listed in the **Well** column is the first well that the sample or control is found in. It is not necessarily the well that failed. Determination of which well failed can be made by examining the results on the NEO Iris.



**Note**: **Well ID** is the control material or the sample ID.

	IRIS	)	Da	tabase Blo	e: 'Galileo_S Immucor R&D dg 7000 Norcross	QL_EN'	<u>.</u>	IRIS
				ABO	DD (QC_ABO	D)		
Plat	e Header Details Plate Name : Al Assay Name : O Printed On : 20 Lot Number : No	BOD QC_ABOD 16-06-21 18:3 ot Defined	1:11		Plate Size : 96 Revision : 2 Printed By : js	8 Wells (8 x 12 inkler	2)	
	Date Saved : 20	16-06-21 18:	31:09		Saved By : jsi	nkler		
Inter	rpretation Warnin	ngs (Plate Va	lid)					
Plea	il Results ise Note - All mod	ified results a	e marked as	RESULT	i			
	Well Type Anti-A Anti-B Anti-AB Anti-A(2) Anti-B(2) Rh-Ctrl Anti-D1 Anti-D1 Diluent Control	Samı 1010 2030: 3011 1023: 2042: 4950: 5170: 5151: 4960:	ble Number 76177261143 33175911096 55174300754 25161081430 66171201583 25190341996 29180961492 19180267187 34170261225	708 613 182 984 092 736 624 115 324	Result Matrix + - - + + +  - + -+ -+ -+ -+		Result Grouping Qualified Qualified Qualified Qualified Qualified Qualified Qualified Qualified	
ndiv Plea	vidual Well Detai se Note - All modi	Is - All Wells ified results ar	e marked as	RESULT				
Well A01 A02 A03 A04 A05 A09 A11 A12 A10	Well ID Anti-A Anti-B Anti-AB Anti-A2) Anti-A2) Rh-Ctrl Anti-O1 Anti-O2 Diluent Control	Result Qualified Qualified Qualified Qualified Qualified Qualified Qualified	A 75 71 11 10 71 11 11 11 11 11 11 11	gglutination \ .500 .700 .100 .900 .900 .700 .400 .100	Alue Well Result + - + + + - - +    	Well Grades 4 - 4 4 4 4 - 4 - 4 - 4 - 4	Pip.Code 1055940505 1055940505 1055940503 1055940505 1055940505 1055940505 1055940505 1055940505 1055940505	
Reag	gent Data							
	Reagent Name corQC Extend Std (1ii Dil dri (10mi) Galileo diluent Anti-DRapid (10mi) H chr (10mi) B-Ceil (10mi) B-Ceil (10mi) Anti-AQ2 (10mi) Anti-AQ2 (10mi) Anti-A (10mi) Anti-A (10mi) Anti-B(2) (10mi)	Lot No. 1224136 496026 5115026 111081 495034 113081 517096 102108 301430 203591 101726 204120	Expiration 2016-07-08 2017-03-25 2018-04-29 2016-06-24 2018-01-29 2016-06-24 2018-01-29 2016-11-20 2017-04-15 2017-03-27 2017-09-23	Barcode 0241901 4960841 0000001 51511911 11117616 4950251 11317610 5170291 1023251 3011051 2030831 1010761 2042661	Well(s) 6136014068811, B1: 7026122532A10, B1: A01, A0: C02, C0 00267187114, 2, B1: 00314199673A09, B0 3034199673A09, B0 3036143262A11, B1: 6108143098A04, C0 7591109661A02, C0 7726114370A01, C0 7726114370A01, C0 77120158309A05, C0	2 2, A03, A04, A05, 3, C04, C05 2, A03, A04, A05, 3, B09, B10, C01, 4 3 2 1 5	A09, A10, A11, A12, B09, I A09, A10 C02, C03, C04, C05	310, B11, B12, CO1,
	Plate Map		a marked as	DESULT				
ull leas	se Note - All modi	fied results ar	e marked as	LOULI				
ull leas	se Note - All modi	fied results ar	e marked as	LOOLI				

#### **Report Naming Conventions and General Features**

In the case of current and sample based reports such as antibody identification and Group and Screen reports respectively, if the assay name selected from the Report Name drop-down list does not match the assay name for which selections are checked in the body of the *Reporting* dialog then a blank grid will be printed on the report.

In the case of plate based reports, if the assay name selected from the Report Name drop down list does not match the assay name for which selections are checked in the body of the *Reporting* dialog then incorrectly formatted data will be printed on the report.

All plates or samples in the *Reporting* dialog are extracted from the filtered list in the main results section and listed in descending chronological order. This allows you to select multiple plates or samples for reporting by clicking the box to the left of the plate/sample ID.

When non-current and current results are selectively combined and printing is attempted, the software displays the *Non-Current Sample(s) Selected* dialog. You must press the **OK** button of the dialog to acknowledge the message.



## Selecting a Large Group of Continuous Entries

To select a large group of continuous entries from the main results screen:

Step	Action
1	Select a plate or sample at one end of the desired group.
2	Press and hold the <b>Shift</b> key on the keyboard.
3	Select the plate or sample at the other end of the desired group.
4	Press the <b>Print</b> button on the software screen.

## Printing Reagents Report from Test Details

To print Reagents reports from Test Details:

Step	Action
1	Press the <b>Details</b> button, for a chosen sample or plate item on the main <i>Results</i> screen under sample or plate view respectively.
	Note: The <i>Reagents</i> report is available using both the sample and plate based results views.
	Select the <i>Reagents</i> tab of the <i>Test Details</i> dialog.
	Note: This version of reagent report lists the reagents used for a given sample or plate. The sample number field displays the sample identification that is currently selected or is defaulted to the top sample in any sample list if a specific sample is not selected, such as in the case of plate based reports. In the case of quality control assays, the sample number field displays the barcoded identification information of one of the reagents undergoing quality control on that quality control plate.
	Test Overview       Results       Resogents       Event Log       Plate Views         Well       Resogent       Lot Number       Qualified       Experision       Last QC       QC Expires       Barcode       Edited       Edited

Step	Action
2	Press the <b>Print</b> button of the <i>Print</i> dialog to print the Reagents report. Press the <b>Cancel</b> button to close the <i>Print</i> dialog. The <b>Number of copies</b> can be selected.
	Print Cancel

## Test Results and Interpretation

For information about Test Results and Interpretation, refer to **Attachment 1 for NEO Iris Operator Manual**.

# **Chapter 9: System Shutdown**

## In This Chapter

This chapter provides a summary of the procedures required to shut down the NEO Iris.

CHAPTER 9:	SYSTEM SHUTDOWN	9-1
Logging Ou	Jt	9-2
Shutting D	own the NEO Iris After Operation	9-3
Extended S	hutdown of the NEO Iris	9-9

## Logging Out

To log out of the NEO Iris:

Step	Action
1	Click the <b>Login</b> button on the Main Menu Bar to log out.
2	The system displays the <i>Login</i> dialog.
	You are now logged out. The screen remains available for the next user to log in.

## Shutting Down the NEO Iris After Operation

#### Purpose

The NEO Iris module and the computer must be completely shut down every day to prevent potential long-term memory leakage and modular drift. Refer to **Chapter 10 – Maintaining the NEO Iris** for information regarding this maintenance task.



**Note**: The maximum duration of the NEO Iris instrument shutdown must be limited to thirty (30) seconds so that the reader lamp and incubators will maintain warmth.

If the NEO Iris remains inactive for a short period (such as overnight), remove all samples, reagents, and unused plates, and store them according to the relevant reagent package inserts.

Because of the washer auto-flush feature, you do not need to perform a full shutdown. The auto-flush feature is the programmed flushing of the washer manifold with 5 ml of system liquid every 30 minutes when the NEO Iris is inactive.



**Note**: If the NEO Iris is in use, auto-flush is not triggered. Therefore, it is important that you monitor the system liquid volume to ensure that sufficient system liquid volume is present, so that the container will not run dry over extended periods of inactivity. Inactivity triggers auto-flush. System liquid monitoring is incorporated into the NEO Iris daily maintenance.



**Note**: You must perform a full system shutdown if you are switching off the NEO Iris instrument power for a prolonged time period (auto-flush will be inactive). Shutdown incorporates flushing of pumps, manifold, and tubing with deionized water to prevent potential salt crystallization build up and blockages. Shutdown also allows you to back up software.

## Summary of Shutdown Strategies

The table below describes the suggested hierarchy of NEO Iris shutdown procedures related to projected time periods of instrument inactivity.

If your goal is to shut down the NEO Iris	You should
to perform the daily maintenance shutdown task.	perform the shutdown procedure, but do not include the Flush System Liquid action (by excluding step 3 of the procedure below).
for a period of time when use of the NEO Iris	perform the shutdown procedure and include the Flush

for testing is not expected in the very near	System Liquid action that leaves deionized water in the
future, but is anticipated to be further in the	tubing. Make a note indicating there is deionized
future (e.g. greater than four days).	water in the tubing to prevent accidental running of
	the instrument before switching the liquid over to PBS.
for a prolonged period (such as taking the	perform the extended shutdown procedure which
NEO Iris out of service and putting it into	includes decontamination and flushing the system with
storage).	air.

#### Procedure

To shut down the NEO Iris:

Step	Action
1	Press the <b>Shutdown</b> button on the Main Menu Bar of the Machine Monitor dialog.
	- SAL
	The shutdown procedure is guided by the Shutdown dialog. Pressing the Shutdown button
	on the Main Menu Bar displays the <i>Shutdown</i> dialog.

Step	Action
2	If the NEO Iris is active then the <b>Shutdown</b> dialog will display with the <b>Ignore</b> button as available, and the running man icon visible. Pressing the <b>Ignore</b> button can have unwanted consequences, as detailed on the dialog. By pressing the <b>Ignore</b> button, the <b>Shutdown</b> and <b>Maintenance</b> buttons will become available.
	Shutdown Procedure         • Please wak until the system is finished If you use [grove to shut down anyhow, processing adely is not ensued.         • If the system is shut down for an extended period of time please fully the spheter and washer with DI water.         • If the system is shut down for an extended period of time please fully the spheter and washer with DI water.         • Please remove all sample racks, reagent racks and plates! If you shutdown with material on board tracking of location changes is not possible.
	Initiate Backup after closing     Shutdown     Operating System
	If the NEO Iris is not active then the <b>Shutdown</b> dialog will display with the <b>Ignore</b> button not available, and the running man icon will not be visible. The <b>Maintenance</b> , <b>Shutdown</b> and <b>Cancel</b> buttons will be available. Note: The <b>Operating System</b> button is only available to Immucor personnel. You can select the <b>Cancel</b> button if you decide to terminate the Shutdown procedure.
	Shutdown Procedure            - Please wat unit the system is frished typou use ignore to statut down androw, processing skety into resulted.             - Do not shut down for an extended period of time please fully the system is shut down for an extended period of time please fully the system is shut down for an extended period of time please fully the system is shut down for an extended period of time please fully out shutdown with material on board.             - Please remove all sample tacks, reagent tacks and plates!         Hypou shutdown with material on board.         Hypou shutdown changes is not possible.
	Initiate Backup after closing     Shutdown     Shutdown     Detroting System     Lancel
3	To run the Flush Liquid System maintenance task: Mote: Do not include this step 3 if performing daily maintenance shutdown.

Step	Action		
	Step	Action	
	1	Select the <b>Maintenance</b> button. This will direct you to the <b>Maintenance</b> screen.	
	2	Select the <b>Flush System Liquid</b> maintenance action from the maintenance list. Refer to <b>Flushing the System Liquid</b> in <b>Chapter 10 – Maintaining the NEO</b> <b>Iris</b> for details of this procedure.	
		Note: Running the Flush System Liquid maintenance task is optional. You should only perform this action if your intention is to leave the shutdown NEO Iris liquid system filled with deionized water. If your intention is to leave the shutdown NEO Iris liquid system filled with PBS, then you should skip the Maintenance button and proceed to remove all sample racks, reagent racks and microplates, and then select the Shutdown button.	
4	To remove	all sample racks, reagent racks and microplates:	
	Step	Action	
	1	Remove all samples from the loading bay and store them according to standard laboratory practices.	
	2	Remove all reagents from the loading bay and store them according to the relevant reagent package inserts.	
	3	Remove all plates from the loading tower and discard all consumed test wells.	
5	To use the	Shutdown button:	
	Step	Action	
	1	Select the <b>Shutdown</b> button. The NEO Iris software completely shuts down.	
Step Action			
--	--		
2 Switch off the NEO Iris module using the power switch on the bottom right side of the NEO Iris module.			

# Extended Shutdown of the NEO Iris

### Purpose

If the NEO Iris is to be taken out of service for an extended period of time, it is recommended that the system be decontaminated, flushed with deionized water, and then primed with air to clear all liquid from the system.



**Note**: If the NEO Iris will be out of service for a short period of time, it is recommended to leave the NEO Iris powered on to allow the auto-flush feature to periodically prime the washer manifold. A sufficient volume of system liquid must be maintained in the system liquid container to prevent the washer manifold from running dry.

## Procedure

To perform extended shutdown:

Step		Action
1	Perform the <b>Iris</b> .	e decontamination procedure as outlined in Chapter 10 - Maintaining the NEO
2	Perform the <b>Iris</b> .	e Flush System Liquid procedure as outlined in Chapter 10 - Maintaining the NEO
3	To prime t	the system with air:
	Step	Action
	1	After completing the Flush System Liquid procedure, disconnect the system liquid tubing connections, including the washer liquid line, from the 20 L system liquid container lid. Rest the ends of the tubing on an absorbent material (e.g. paper towel).
	2	Run the <b>ClearTube</b> assay. This will prime the system with air to remove all liquid in preparation for storage.
	3	Once the assay is complete, remove both system liquid containers and discard the deionized water. Return the empty containers to their proper position on the cabinet drawer.

	4	Reconnect the system liquid and washer tubing to the empty 20 L system liquid container lid.
4	Enter the S Iris software centrifuge r	hutdown dialog as described in this chapter. Select the <b>Shutdown</b> button. The NEO e completely shuts down. Switch off both the NEO Iris instrument and the module using the on/off switches for both systems.
		<u>Note</u> : When bringing the NEO Iris back into service following an extended shutdown, perform the decontamination procedure as outlined in <b>Chapter 10</b> - <b>Maintaining the NEO Iris</b> .

# **Chapter 10: Maintaining the NEO Iris**

# In This Chapter

This chapter describes the detailed procedures for maintaining the NEO Iris.

CHAPTER 10: MAINTAINING THE NEO IRIS	10-1
Maintenance Overview	10-2
Daily Maintenance Tasks	10-6
Weekly Maintenance Tasks	10-28
Monthly Maintenance Tasks	10-46
As Needed Maintenance Tasks	10-61

# Maintenance Overview

### Maintenance and Verification Action Status Screen

The *Maintenance and Verification Action Status* screen lists all of the maintenance actions that are available on the NEO Iris, along with the current state, the required interval, and the date the next action must be performed. You can highlight and perform multiple actions at one time.

To access the Maintenance and Verification Action Status screen, click the Maintenance button.



The **Maintenance and Verification Action Status** screen lists all of the software driven maintenance actions that are available on the NEO Iris, along with the current state, the required interval, and the date the next action must be performed. You can highlight and perform multiple actions at one time.

Action	State	Interval Next Action 🔺	
			-
ABO Reagent QC	Due	1 day(s) 2009-08-15	
Service Function 1	Due	0 day(s) 2009-08-18	
Clean Instrument	Due	1 day(s) 2009-08-18	_
Check Pipettor Reference	Due	1 day(s) 2009-08-18	Ξ
Verify Washer	Due	1 day(s) 2009-08-18	
Verify Reader	Due	1 day(s) 2009-08-19	
Pipettor Self Check	Due	1 day(s) 2009-08-19	Ŧ
Flush Liquid System	Due	0 day(s) 2009-08-19	-
Pipettor Verification Test	Due	0 day(s) 2009-08-19	<b>_</b>

The NEO Iris software defines monthly as 30 days. You can successfully perform any maintenance action **before** the required time interval date. If you do, the next required time interval date for that action remains automatically calculated from that last action date. For example, if you performed a 7-day action after only 4 days, the next required date is calculated as 7 days from this latest action date. **Daily** is defined as a 24-hour period, such as 9 AM to 9 AM.

The system lists the state of maintenance actions as **OK**, **Due**, or **Pending**. **OK** means no action is required. The next action for that item is designated for the future. **Due** means that this action must be taken now and the required date is overdue. The text for a **Due** item is in **red** font. Actions with an interval of 0 days do not prevent assays from being run, and it is not necessary to perform these actions

when their state is **Due**. **Pending** means that this action is in process and awaiting completion. Once an action is completed, the status changes from **Due** to **OK**. A status of **Not Passed** can be associated with a task if that given task failed. An example of such a task is **Pipettor Self Check**.

To run a particular maintenance action, you can select (highlight) the given action and then press the **Start** button on the **Maintenance and Verification Action Status** screen. The on-screen instructions will then be displayed for that action. Press the **Done** button to exit the **Maintenance and Verification Action Status** screen.



Press the **Continue** button of the **Maintenance and Verification Preparation** screen to begin the maintenance action, or press the **Cancel** button to abort the action.

You can view an archive list of a given software driven maintenance action by selecting the required action on the **Maintenance and Verification Action Status** screen and then selecting the **View** button. The **Maintenance and Verification Archive** screen is then displayed for that required action which lists each chronological event.

Action Pipettor Sell Direck	Maintenance and Ver	lication Archiv	re .	
Anay	Date +	State	Operator	I
PDwok	2010/02/25 13:56:06	OK .	Donna	
PDeck	2010-02-24, 13:46:47	ÓK.	Yolarida	1
PDreck	2010/02/24, 11:39:56	OK.	Yolanda	1
PCheck	2010/02/24, 11:37:20	OK.	Yolarda	-
POwek.	2010/02-24, 10:54 29	OK.	Volanda	Ξ
POwek	2010/02/23, 14:11:15	0K	Tanvor	
PDreck	2010-02-23, 07-57-03	OK.	Stacy	
PDwck	2016/02/23, 07:51:08	OK.	Staty	3
PCheck.	2010/02/22, 13:36:36	OK.	Tievor	
PCheck.	2010-02-22, 13:30:03	OK.	Tievor	×
ок			Delete	Detail

You can view the event details by selecting a particular event from the **Maintenance and Verification Archive** screen and then pressing the **Detail** button. The event detail screen is then displayed. You can delete an archived log entry by highlighting a particular entry on the archive screen and pressing the **Delete** button, but you must have the necessary security level to be able to delete.



Various aspects of an event detail, other than **General Information**, can be viewed by selecting the other associated tabs (Assay Protocol, Plate Events, Sample ID, Raw Results and Flags).

Press the **OK** button of the event detail screen to return to the **Maintenance and Verification Archive** screen. Press the **OK** button of the **Maintenance and Verification Archive** screen to return to the **Maintenance and Verification Action Status** screen.



**Limitation:** Maintenance actions for the NEO Iris verify that specific modules of the instrument are functioning at the required specifications. The actions described are critical to the NEO Iris assay performance.

If a maintenance item is not successfully performed for a module within the required time interval, then any assays that use that module will not run. For example, all assays use the reader. If you do not perform the daily **Pipettor Self Check** task then all assays are locked down and no assays can run.



**Limitation:** Before using any cleaning or decontamination methods except those recommended by the manufacturer, users should check with the manufacturer that the proposed method will not damage the equipment.



**Warning:** If a maintenance action requires the manual entry of a plate ID, **do not** input an ID with 2 digits of alpha characters followed by numeric characters. Because the software may recognize this ID scheme as that of a sample testing plate, entry of a plate expiration date will be required. Entry of a plate expiration date is not part of the instructions for maintenance actions.

You can log documentation of maintenance actions onto the four maintenance forms: NEO Iris Daily Maintenance Record, NEO Iris Weekly Maintenance Record, NEO Iris Monthly Maintenance Record, and PipTest Chart. Master copies of these four record sheets are located in **Appendix B – Maintenance Records**.

# Daily Maintenance Tasks

# **Daily Tasks**

The following tasks are to be completed each day:

- Perform a NEO Iris module and PC shutdown
- Filling the system liquid
- Emptying the common waste container
- Clean the instrument
- Check Pipettor Reference
- Complete Reagent Quality Control (QC)
- Complete QC3\_Cell Quality Control (QC)
- Pipettor Self Check

## PC and NEO Iris Module Shutdown

### Purpose

The purpose of shutting down, and immediately powering up, the NEO Iris module and PC is to reset all volatile module memory and prevent potential long-term computer memory leakage.

### Procedure

To shut down and power up the NEO Iris module and PC:

Step	Action
1	Shut down the PC according to the daily routine procedure described in <b>Chapter 9 – System Shutdown</b> .
2	Power off the NEO Iris module.
3	After step 2 above is completed, <u>immediately</u> power on the NEO Iris module and PC according to the procedures described in the starting up section of <b>Chapter 5</b> – <b>Instrument Start-Up</b> .

### Automatic Deletion of Outdated Plates

At the beginning of full initialization, the following outdated information is automatically deleted from the database:

- plates from the plate list, greater than forty eight (48) hours old
- plates from maintenance actions with greater than two (2) plates per assay (leaving 2 plates in the database after the automatic deletion), and
- event log entries, greater than ninety (90) days old.

As a result, the manual deletion of historic sample and maintenance plate entries is not required on the part of the operator.

## Filling the System Liquid

Both of the system liquid containers are mounted inside the cabinet on a sliding drawer that allows easy access to the liquid containers and connections when pulled forward. This section describes how to check the 10-liter refill system liquid container and fill it with system liquid.

### Procedure

To fill the 10-liter System Liquid Container:

Step	Action
1	Open the cabinet doors and pull out the drawer to access both of the system liquid containers.
2	Check the system fluid level of the 20-liter system liquid container. If required, follow the procedure below to refill the container with system liquid.
3	Remove the cap of the 10-liter refill system liquid container and fill it with system liquid. A funnel can be used. Alternatively, the 10-liter refill container can be disconnected and filled off-line. Reconnect the 10-liter refill container after filling this container off-line. Allow approximately five (5) minutes for gravitational equilibration of liquid levels into the 20-liter container.
	Left: 20-liter system liquid container. Right: 10-liter refill container.

Step	Action
4	Replace the cap of the 10-liter refill system liquid container, to prevent contamination.
	<b>Warning:</b> Results can be adversely affected if the system liquid container is filled with anything other than PBS (e.g. de-ionized water).
5	Click the Washer button on the Machine Monitor. The system displays the Wash Buffers dialog.
	The <b>Buffer ID</b> name will automatically be displayed as PBS.
6	The green Fill Level check mark (✓) appears in the pump area of the Wash Buffers dialog when the system liquid container has adequate fluid. This must be the case in order to prime the pump.         Image: Note: The pump will automatically prime itself if Initialization is run.         If needed, press the Initial Prime button. This primes the pump with PBS. The In use circle turns red.         Image: Note: If the 20-liter container is filled, but the Initial Prime button is not available, click in the Buffer ID field for the pump. This activates the Initial Prime button
7	Click Done to leave the dialog
1	Click Duile to leave the ulaiby.

Step	Action
8	When you have completed steps 1–7, ensure that the sliding drawer is fully pushed back into
	the cabinet and close the doors.
	<b>Warning:</b> Filling the container with the incorrect fluid adversely affects Solid Phase
	Red Cell Adherence reactions.



**Limitation:** Immucor requires the use of phosphate buffered (approximately 15mM) isotonic saline, pH 6.5-7.5 (PBS), on the NEO Iris system. Reactions between an antibody and its antigen may be weakened if acidic or unbuffered saline is used.

Using saline and/or deionized water in PBS preparation from sources with systems in place to control proliferation of microbes helps to reduce the chance for microbial bioburden on the system. Excessive microbial bioburden can cause degradation of system or assay performance.

### Use of Cubed Saline

As an alternative method of filling, a commercially obtained cube of unbuffered saline (that must be buffered prior to its use as per limitation above ) can be directly connected (using the supplied in-line connector) to the 20-liter container, instead of using the 10-liter refill container. In this instance, the commercially obtained cube of PBS must be elevated on the drawer to at least the same level as the 10-liter refill container so that passive gravitational draining equilibration can take place. The cube must be located in the space on the drawer to the left of the common waste container and requires an extended length of tubing to connect the cube to the detachable in-line connector.

## **Emptying the Common Waste Container**

The purpose of emptying the common waste container is to maintain empty space in the waste container for the collection of waste fluids (into the waste container).

You must perform this task at least every day or when the common waste container is full. The waste container status indicator, located on the bottom of the screen, monitors the liquid level in the container. It is normally displayed with a yellow background. When the waste container is full, the status indicator is displayed with red background.



Waste container status indicator in a capacity-for-waste condition

### Required Tool(s)

The following tool(s) are required to complete this task:

• Waste shuttle container

### Procedure

To empty the common waste container:

Step	Action
1	Obtain the empty waste shuttle container before starting this procedure.
2	Open the cabinet doors and pull out the drawer to access the common waste container.

Step	Action
3	Disconnect the female quick-connector of the draining tube from the lid of the common waste container.
4	Connect the female quick-connector at the end of the draining tube to the male quick-connector on the waste shuttle container cap. You will hear a "click" when the connection is made.
	Female quick-connector connected to the waste shuttle container cap
	The liquid waste drains into the waste shuttle container.

Step	Action
5	When the waste transfer is complete or the shuttle container is full, disconnect the draining tube from the waste shuttle container cap.
6	Reconnect the draining tube to the male quick-connector on the lid of the common waste container. You will hear a "click" when the connection is made.
7	Push the drawer back into the cabinet and close the doors.
8	Dispose of the waste liquid following your standard laboratory practice. <u>Warning</u> : Blood samples, liquid waste, used micro-well strips, and consumed liquid reagent vials contain potentially bio-hazardous material.
	Always wear protective gloves and clothing when handling blood samples, liquid waste, micro-well strips, or liquid reagent vials. All blood samples, liquid waste, micro-well strips, and liquid reagent vials must be discarded following the standard practice of the laboratory.
	All blood products must be treated as potentially infectious. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

# **Cleaning the Instrument**

### Purpose

The purpose of cleaning the instrument is to remove biohazardous contamination on the NEO Iris caused by spills or debris from reagents or samples.

NEO Iris initialization is included as part of the cleaning process.

### Required Tool(s)

The following tool(s) are required to complete this task:

• Alcohol wipes or 0.25% sodium hypochlorite (NaClO) as the recommended cleaning solution *Refer to Monthly Maintenance- Decontaminate Tubings for the specifications and preparation of 0.25% NaClO* 

### Procedure

To clean the instrument:

Step	Action
1	Press the Maintenance button on the Main Menu Bar.
2	Select <b>Clean Instrument</b> from the action list.
3	Press the <b>Start</b> button.

4	Follow the on-screen instructions.
	Please perform the following preparation steps:
	Clean Instrument <housing and="" cover="" front="">:</housing>
	1. Clean insrument housing and front cover using a paper towel moistened with 70% isopropanol or recommended cleaning solution.
	2. Repeat cleaning using deionized water.
	Clean Instrument <plate loading="" tower="">:</plate>
	1. Clean the Loading Tower using a paper towel moistened with 70% isopropanol or recommended cleaning solution.
	2. Repeat cleaning using deionized water.
	Clean Instrument <loading bays="">:</loading>
	1. Clean the sample and reagent loading bay using a paper towel moistened with 70% isopropanol or recommended cleaning solution.
	2. Repeat cleaning using deionized water.
	Clean Instrument <pipetting positions="">:</pipetting>
	1. Turn off the instrument
5	Press the <b>Continue</b> button below the instruction screen after the NEO Iris initialization is completed.
	·

# **Checking Pipettor Reference**

### Purpose

The purpose of checking the pipettor reference is to

- Verify that pipettor alignment is correct prior to running assays.
- Correct pipettor alignment after a probe has been bent. You should first straighten the probe.
- Confirm that alignment is acceptable after probe exchange.

### Procedure

To check pipettor reference:

Action
Press the <b>Maintenance</b> button on the Main Menu Bar.
Select Check Pipettor Reference from the action list.
Press the <b>Start</b> button.
Follow the on-screen instructions.
Please perform the following preparation steps:
<ul> <li>Check Pipettor Reference <probe alignment="">:         <ol> <li>Be sure the pipettors are idle prior to checking the reference positions.</li> <li>If the probes need slight adjustments, gently bend to the target.</li> </ol> </probe></li> <li>Marning: This maintenance task can only be performed when the NEO Iris is not active. You should only use this function for minor corrections. If a probe is badly bent, you should straighten it or replace it. Resetting a badly bent probe can result in crashes as the hardware tries to move to positions outside of the mechanical boundaries of its reach.</li> </ul> <li>Marning: You can open the instrument hood during this maintenance action to access the probes. When accessing the instrument deck, you must not move the device of the instrument deck.</li>

Step	Action
5	Press the <b>Check Reference</b> button.
6	Select the Reference Position button. The instrument defaults to checking the left arm first.         Image: Check Pipettor Reference Position         Image: Pipettor arm Image: Position         Image: Provide the provided of the problem of t
7	If corrections are required for the left arm in the X (left-right) or Y (front-back) direction, gently bend the probe(s) to the target.
8	All probes must be centered approximately 0.75 mm above the target surface and all must be horizontally even. The probes must not touch the target. If up-down adjustments are necessary, you must contact Technical Support to determine further action.
9	After the probes are aligned correctly, click the <b>Position OK</b> button. You can now check the right arm.
10	If there was a problem aligning the left arm, you should click the <b>Position Wrong</b> button. By clicking the <b>Position Wrong</b> button, the process will be returned to step 6 above.
11	Select the <b>Right</b> arm and click <b>Reference Position</b> . The right arm moves to the reference position target and the probe will move to the previous Z-position.           Mote:         Corrections are not necessary if the probe tip is within the inner circle of the reference target.
12	If corrections are required for the right arm in the X (left-right) or Y (front-back) direction, gently bend the probe to the target.

Step	Action
13	Inspect the probe. The probe must be centered approximately 0.75 mm above the target surface. The probe must not touch the target.
	If up-down adjustments are necessary, you must contact Technical Support to determine further action.
14	After you have correctly aligned the probe, click the <b>Position OK</b> button.
15	If there was a problem aligning the right arm, you should click the <b>Position Wrong</b> button. By clicking the <b>Position Wrong</b> button, the process will be returned to step 11 above.
16	The <b>Close</b> button should now be selected. This will initialize both arms and return you to the instruction screen.
	Continue should be selected to complete the task and return to the Maintenance screen.

# **Completing Reagent Quality Control (QC)**

### Purpose

The purpose of completing the reagent QC is to:

- Verify that the agglutination reagents are reacting appropriately.
- Verify that the centrifuge shake sequence is performed correctly.

You must perform this action every day, or prior to the use of all agglutination reagents.

If the agglutination reagents are not used for one or more days, then it is not necessary to perform the reagent QC until just prior to when those reagents are required for blood sample testing. In this instance, the reagent QC maintenance action will be listed on the Maintenance and Verification Action Status screen as **Due** until the QC is performed.



**Warning:** Immucor requires that the quality control (QC) assay for agglutination tests is performed on the system every day, or prior to the use of agglutination reagents. The primary objective of QC is to ensure that the loaded reagents are satisfactory.



**Warning:** Repeated QC failures in absence of a reagent related cause can be an indication of conditions that have impacted instrument system performance (i.e. pipetting, centrifugation, incubation, etc), or that an incorrect system liquid (deionized water rather than PBS) may have been loaded.



**Note:** The procedure described in this section is provided as a standard method which fits with most of the local QC requirements. However, if this Reagents QC is not appropriate for your region, please contact your local Immucor Representative to accommodate the QC procedure to your local requirements and regulation needs.

### Reagents

Refer to the NEO Iris **Assay Reagent Component Grid** of Attachment 1 for NEO Iris Operator Manual for the list of reagents used with reagent QC.

The reagent QC assay includes agglutination reagents for ABO- and Rhesus D-typing and for Rh- and Kell- phenotyping.

### Procedure

To complete the reagent QC check:

Step	Action
1	Load an unused ABO plate by the methods described in <b>Chapter 6 – Instrument Testing Operation</b> .
2	Load all of the necessary reagents for the <b>QC</b> assay, as listed in <b>Attachment 1 for NEO Iris</b> <b>Operator Manual</b> , by the methods described in <b>Chapter 6 – Instrument Testing Operation</b> .
3	Press the Maintenance button on the Main Menu Bar.
4	Select the <b>Reagent QC</b> item on the <b>Maintenance and Verification Action Status</b> screen, and then press the <b>Start</b> button. This will take you to the <b>Maintenance and Verification</b> <b>Preparation</b> instructions for <b>Reagent QC</b> .
5	Press the <b>Continue</b> button on the <b>Maintenance and Verification Preparation</b> screen to access the <b>Resource Overview</b> window.
6	Select the <b>QC</b> line on the <b>Resource Overview</b> window to determine if any additional resources are required. Load any additional resources required.
7	Press the <b>Start</b> button of the <b>Resource Overview</b> window to begin the <b>QC</b> assay when all resource requirements are satisfied.
8	Once <b>QC</b> is complete, remove the plate from the Plate Loading Tower and discard it.

# Completing QC3\_Cell Quality Control (QC)

### Purpose

The purpose of completing QC3\_Cell is to:

- Verify that the 3\_Cell reagents are reacting appropriately.
- Verify that centrifugation of 3\_Cell assays is performed correctly.

#### **Required Interval**

You must perform this action every day, or prior to the use of 3\_Cell reagents.

If 3\_Cell reagents are not used for one or more days, then it is not necessary to perform the QC3\_Cell QC until just prior to when those reagents are required for blood sample testing. In this instance, the **3\_Cell Plate QC** maintenance action will be listed on the **Maintenance and Verification Action Status** screen as **Due** until the QC is performed.



**Note:** Immucor provides a QC for the 3\_Cell reagents only. Other Capture screening methods (using 2\_Cell, 4\_Cell or Pool\_Cell reagents) always use an internal Negative control and an internal Positive control (containing a mixture of antibodies) with each run. Therefore an additional daily QC is not required for those as it is built into the test system. Please contact your local Immucor Representative if you need more details and/or if you want to discuss the possibility of using other external QC material on a daily basis with the 2\_Cell, 4\_Cell or Pool\_Cell screen assays.

#### Reagents

Refer to the NEO Iris **Assay Reagent Component Grid** of Attachment 1 for NEO Iris Operator Manual for the list of reagents used with QC3\_Cell.

#### Procedure

To complete the QC3\_Cell quality control:

Step

Action

Step	Action
1	Load an unused Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (3) strip into a Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (3) white frame at strip position one (1). Load the prepared plate by the methods described in <b>Chapter 6 – Instrument Testing Operation</b> .
	Warning: For the QC3_Cell assay, the correct position for the Capture-R® Ready- Screen® (3) strip in a Capture-R® Ready-Screen® (3) white plate frame is at strip position one (1). The instrument will not check for the correct placement of the Capture-R® Ready-Screen® (3) strip in the white plate frame. You must make sure of the correct placement of this strip yourself. Incorrect placement of this single strip can result in test material being pipetted onto the pipetting station surface, an overflow error in the washer, and a failed quality control. Unused Capture-R® Ready-Screen® (3) strips should not be present in any of the other positions of the plate. It is not required that strip positions two (2) to twelve (12) of the white plate frame be occupied by any strips at all. Those positions can be left empty, but all strip positions should be activated in the software. The Strip Selection tab in the Plate Loading Tower dialog must indicate that all twelve (12) strips are available even though only the strip at position (1) is present.
2	Load all of the necessary reagents for the QC3_Cell assay, as listed in Attachment 1 for NEO Iris Operator Manual, by the methods described in Chapter 6 – Instrument Testing Operation.
3	Press the Maintenance button on the Main Menu Bar.
4	Select the <b>3_Cell Plate QC</b> item on the <b>Maintenance and Verification Action Status</b> screen, and then press the <b>Start</b> button. This will take you to the instructions for <b>3_Cell Plate QC</b> .
5	Press the <b>Continue</b> button on the <b>Maintenance and Verification Preparation</b> screen to access the <b>Resource Overview</b> window.
6	Select the <b>QC3_Cell</b> line on the <b>Resource Overview</b> window to determine if any additional resources are required. Load any additional resources required.
7	Press the <b>Start</b> button of the <b>Resource Overview</b> window to begin the <b>QC3_Cell</b> assay when all resource requirements are satisfied.
8	Once <b>QC3_Cell</b> is complete, discard the used Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (3) strip, but save the white frame for future use.

# Performing the Pipettor Self Check

### Purpose

The purpose of performing the pipettor self check is to verify that the probe wash tower rinse pumps are functioning correctly and that the probes are aspirating and dispensing fluid correctly. Both the four probe arm and the single probe arm are verified.

### Required Tool(s)

The following tool(s) are required to complete this task:

• Pipettor self check tool

### Procedure

To perform the pipettor self check:

Step	Action
1	Select the Pipettor Self Check action from the Maintenance and Verification Action Status
	list and then press the <b>Start</b> button.



Step	Action
4	Close the door of the Plate Loading Tower to initiate an automatic reading of the tool
	barcode.
	Note: The tool will be taken to the plate barcode scanner by the transport so that the barcode can be read. The tool will then be returned to its assigned position in the Plate Loading Tower. Wait until the tool has been returned to the Plate Loading Tower before proceeding to step <b>5</b> .
	Barcode
5	Press the <b>Continue</b> button on the <i>Pipettor Self Check</i> instruction screen.
	Note: Allow step 4 to complete before performing this step 5.
6	Select the <b>PCheck</b> assay on the <i>Resource Overview</i> window.
7	Press the <b>Start</b> button on the <i>Resource Overview</i> window to begin the PCheck assay.
	The software analyzes the results and sets a pass or fail status. Refer to
	Chapter 11 – Troubleshooting the NEO Iris for details regarding information about the fail
	status.
	Note: No Liquid/Not enough liquid flags are generated and appear in the Log List during processing. This is expected and necessary.
8	After completion of the PCheck assay, remove the pipettor self check tool from the plate
	loading tower and then also from the plate carrier.

Step	Action
9	Rinse the cavities only of pipettor self check tool with tap water and then immediately dry all of its surfaces carefully with a paper towel.           Note:         Make sure not to damage the barcode affixed to the tool.
10	Store the pipettor self check tool in a dry secure place.

# Weekly Maintenance Tasks

# Weekly Tasks

The following tasks are to be completed weekly:

- Cleaning the Reader Mirror and Renewing Flatfield Images
- Creating Archives
- Cleaning the Pipettor Wash Towers
- Cleaning the Common Waste Container
- Inspecting the Syringes

# **Cleaning the Reader Mirror and Renewing Flatfield Images**

### Purpose

The purpose of cleaning the mirror and renewing flatfield images is to:

- To remove spills or dirt on the mirror hampering the reader functionality.
- Take background images to eliminate dust and debris from interfering with the interpretations.

### Required Tool(s)

The following tool(s) are required to complete this task:

- Alcohol wipes
- Clean dry absorbent laboratory disposable wipes
- Flashlight

### Procedure

To clean the reader mirror:

Step	Action
1	Power off the NEO Iris module.
2	Refer to the instructions provided to clean the mirror in the <b>As Needed Maintenance Tasks</b> section of this chapter.

To renew flatfield images:

Step	Action
1	Power on the NEO Iris module.
	<b>Limitation:</b> The instrument must be switched on for at least thirty (30) minutes prior to take flatfield images so that the reader lamp can warm up.
2	Press the <b>Maintenance</b> button on the Main Menu Bar.

Step	Action
3	Select Take Flatfield Images from the action list.
4	Press the <b>Start</b> button.
5	Follow the on-screen instructions for taking the Flatfield Image.
	Maintenance and Verification Preparation Please perform the following preparation steps:    Please perform the following preparation steps:   Please perform the reader is idle. Select Take Elablish
	2. Check the flatfield images for spills that may have occurred.
	3. If a large spill is present, remove the cover plate of the reader module and clean the mirror.
	4. Retake the flatfield image. If the spill is still visible remove the panel to the reader module and clean the light box at the top of the module.
	5. Retake the flatfield image and if the spill has been cleaned continue to the next step.
	Warning: You should not initiate this maintenance         function if an assay is scheduled to use the reader.         Leaving a microplate in the reader during this action         affects the analysis of results.
6	Press the <b>Take flatfield</b> button after confirming that a plate is not in the reader.
	Take flatfield
	The system takes one flatfield image for "Default Setting" analysis used on the NEO Iris. While
	the system takes the image, look for any large spills that may have occurred on the mirror or
	the light diffuser.

Step	Action
7	Once the system has finished taking the flatfield image, the system makes the <b>Done</b> button available in the <i>Reader Flatfield Image</i> window.
	Press the <b>Done</b> button.
	If there are any significant spills, the mirror and the light diffuser need to be cleaned. Refer         to the instructions provided for these tasks in the As Needed Maintenance Tasks section of
	this chapter.
8	Repeat steps 5 and 6 to verify that the spill is completely cleaned.

### **Creating Archives**

### Purpose

The purpose of creating archives is to reduce the physical size of the database and to preserve historical results in an electronic format. Additionally, during the archive process, excess space created by the process of deleting and modifying data is removed. Performing this maintenance task can minimize software errors caused by database corruption and data fragmentation.

### The Incorporated Warning System

A yellow warning triangle will be displayed on the **Utilities** button of the Main Menu Bar if the configured required time interval for the scheduled next archive of results has been exceeded.



Pressing the **Utilities** button (with the warning displayed) gives access to the **Utilities** window. A second yellow warning triangle is subsequently associated with the **Archive** tab of the **Utilities** window.

Event log	🔔 Archive	Statistics	Reports				
-Drive space: Final archive location							

The yellow warning triangle is removed after the required archiving is completed which then subsequently restarts the configured required time interval. Refer to **Chapter 4 – Security** for details regarding **Archive Configuration**.



**Warning:** Failure to archive results according to the configured time schedule interval, at a minimum, can result in the accumulation of excessive result data on the instrument computer which can lead to the slowing down of the speed of the computer processing and activity.



**Warning:** You must not enter the *Archive* tab of the *Utilities* window while the NEO Iris is processing, otherwise there is a risk that conversion of results could be interrupted and that results could need to be re-sent.

### Procedure

To create archives:

Step	Action				
1	Make sure that the selected backup medium is attached to the PC (i.e. external USB drive) or insert a DVD+R disk into the disk drive. If DVD disk is used, wait until the LED on the drive is no longer illuminated before moving to step 2.         Image: Note: Do not use CD or DVD-R disks.         Image: Note: Only one archive is allowed per disk.				
2	Select the <b>Utilities</b> dialog from the Main Menu Bar.				
3	Select the Archive tab of the Utilities window.				
4	Verify that the space estimate required for this archive is less than the Free space in both the Final archive location and the Local archive location.				
Event log Archive Statist	ics Reports			OK	
---	-----------------------	------------------------	--------------------------	------------------	----------------------
Drive space: Final archive I Capacity: 702.82 MB	Current	action :		0%	
Used: 0 Bytes		Total :		0%	
Drive space: Local archive Capacity: 35.62 GB	ocation				
Used: 5.45 GB					
Estimates					
Space required: 259.1 MB Status					
Oldest existing data: 2009 Last archive date: 2009	06-02 10-09				
Next archive date: 2009 Working database: Galile	10-16				
Action Archive -	nd images				
C Copy	Archive logs			Go!	
10, 100,550					
	Note: Before p	performing the	archive process	, the instrum	ent calculates the
	amount of spa	ace required to	archive the dat	a and the am	nount of space
	available on th	ne archive med	lia and on the lo	ocal drive. If t	here is less space
	the archive me	edia or on the	local drive than	is needed fo	r the archive the
	instrument wil	l alert vou with	he the following of	lialog. The di	ialog will indicate
				lialog. The u	
	the drive that	nas insufficien	t space is the fin	ial archive m	edia or the local
	temporary fold	der. In this inst	ance, the archive	e will not be	performed. Clear
	dialog by pres	sing the <b>OK</b> b	utton and verify	the proper n	nedia is in the driv
	the message p	persists, contac	t Immucor Tech	nical Support	ī.
Warning: Availa	ble Space				
There may	not be enough free sp	ace on the Final targe	t media for this Archive		
Estimated	Required: 377.24 MB				
Available:	-				
Free some	space before starting	the Archive			
		OK			
		~~			
	Note: The soft	tware will auto	matically name t	the archive a	ccording to the d
	range to be ar	rchived. The yy	yymmdd_YYYYM	IMDD format	t is used where
	yyyymmdd is t	the date corre	sponding to the	oldest data k	peing archived an
	YYYYMMDD is	the date corr	esponding to the	e most recen	t data being archi
	according to t	he number of	dave kent in the	working dat	abase (i.e. current
	data must		aurod to ke he	t) In sector	
	uate – numbe	r of days confi	gurea to be kep	i). In cases w	inere more than c
	archive is crea	ted with the sa	ame date range	a "_N" will be	e appended to

6	The <b>Progress</b> of the archive is displayed as it proceeds. Once complete, the scroll bar can be used to go back and view the full sequence.
	<b>Note:</b> Do not attempt to stop the archive or turn off the PC while the archive is in progress. This may result in data corruption or data loss.
	Progress Current action : 56% Total : 24%
	Moving image file D:(Galleo)images/RVREIFY_488500:1_0.bbc/         Moving image file D:(Galleo)images/RVREIF/488500:1_0.bbc/         Moving image file D:(Galleo)images/RVREIF/488500:0_0.bbc/         Moving image file D:(Galleo)images/RVREIF/485007_0.bbc/         Moving image file D:(Galleo)images/RVREIF/485007_0.bbc/         Moving image file D:(Galleo)images/RVREIF/485007_0.bbc/         Moving image file D:(Galleo)images/RVREIF/485007_0.bbc/         Moving image file D:(Galleo)images/RVREIS0300_4851174_0.bbc/         Moving image file D:(Galleo)images/RVREIS03000_4851174_0.bbc/         Moving image file D:(Galleo)images/RVREIS03000_4851164_0.bbc/         Moving image file D:(Galleo)images/RVREIS0300_4851164_0.bbc/         Moving image file D:(Galleo)images/RVREIS0300_4851164_0.bbc/         Moving image file D:(Galleo)images/RV
7	When the archive process is completed, a confirmation dialog is displayed stating <b>Archive complete!</b> Press the <b>OK</b> button of the <b>Archive complete!</b> message. Do not remove the archive media from the drive at this time.
	Archive complete!
8	Press the <b>OK</b> button to close the Utilities dialog.

If a DVD was used, remove the archive disk from the PC disk drive and store it in a safe place, according to the recommended storage conditions of the disk manufacturer.



**Note:** It is recommended that the disk is removed only when the disk drive LED is no longer illuminated following the display of the message. This ensures that the process has gone to completion.

### **Rewriting Archives**

9

The rewriting of archives to a disk is triggered if there is a problem with the burning of an archive to a disk and the one (1) permitted retry attempt fails. A backup copy of the archive is stored on the harddrive of the PC in such cases.

The re-write process has no limits on the number of times that it can be executed.

To rewrite archives:

Step	Action		
1	Make sure that a DVD+R disk is in the disk drive. Wait until the LED on the drive is no longe illuminated before moving to step 2.		
	Note: Do not use CD or DVD-R disks.		
	Note: Only one archive is allowed per disk.		

Step	Action		
2	If any archives already exist on the hard drive, then when you select the <i>Archive</i> tab, the <b>Unburnt archives located!</b> dialog is displayed.		
	Unburnt archives located!		
	Re-write archive Re-write later		
	Press the <b>Re-write archive</b> button to display a list of all backed up archives on the <b>Select</b> archive to retry dialog.		
	Note:If the Re-write later button is pressed, the Unburnt archives located!dialog is closed and the Archive dialog is then displayed. In this case, a Re- write Archive button will be displayed on the Archive dialog where you can		
	attempt the re-write at anytime.		
	U Archive  So! Go!		
	When all harddrive stored backups have been copied to disk sucessfully, the <b>Re-write</b> <b>Archive</b> button will no longer be displayed on the <b>Archive</b> dialog.		
	Go!		

Step	Action
3	Select the archive to be written to disk from the list and then press the <b>OK</b> button.
	Name     Comments/Notes       Physical archives
	20090929_20100222-3
	OK Cancel
	The <b>Archive re-writer</b> dialog is displayed after the <b>OK</b> button is pressed.
	Alternatively, press the <b>Cancel</b> button to cancel the archive rewrite process.

Step	Action
4	Press the <b>Write Disc</b> button on the <b>Archive re-writer</b> dialog to start the writing of the selected harddrive archive to the disk. Alternatively, press the <b>Close</b> button to close the <b>Archive re-writer</b> dialog.
	Archive re-weiker [NMSDVD Burning SDK for .NET - Version 1.004]     IXI       Device     EHLDT-ST DVD-RAM GSA-H60L RS0C       Image: Strategy school (Strategy school
	After pressing the Write Disc button and the process subsequently completes, a message is displayed indicating that the process is complete and the selected backup copy on the hard
	drive is removed.       Sanguin. Ldr. DataWriter       Image: Write and Verify Complete!
	Press the <b>OK</b> button of the <i>Write and Verify Complete!</i> message.
5	Remove the archive disk from the PC disk drive and store it in a safe place, according to the recommended storage conditions of the disk manufacturer.
	Note: It is recommended that the disk is removed only when the disk drive LED is no longer illuminated following the display of the completion message. This ensures that the process has gone to completion.

## **Cleaning the Pipettor Wash Towers**

### Purpose

The purpose of cleaning the pipettor wash towers is to reduce tower contamination.



Note: The NEO Iris must be switched OFF.

### Required Tool(s)

The following tool(s) are required to complete this task:

• 0.25% sodium hypochlorite (NaClO) as the recommended cleaning solution

Refer to Monthly Maintenance- Decontaminate Tubings for the specifications and preparation of 0.25% NaClO

### Procedure

To clean the pipettor wash towers:

Step	Action		
1	Prepare a concentration of 0.25% NaClO working solution.		
	Note: Approximately 50 to 60 mL of solution is required for this weekly task. This task can be performed at the same time as the weekly decontamination of the waste container and the monthly decontaminate tubings task, if convenient.		
2	Manually move the pipettor arms away from both probe wash towers and reposition them over a loading bay by either moving both arms over the 14-lane Bay, or by moving the left arm over the 14-lane Bay and the right arm over the 5-lane Bay.		
3	Dispense the working solution of the recommended cleaning solution into the probe tip rinse cup of both towers and allow it to drain, or be added directly into the surrounding waste cup until filled. The solution should not be added to the overflow trough and tubing. The volume of solution required to completely fill each tower can range between 25 to 30 mL (approximately 50 to 60 mL in total).		

4 Allow the pipettor wash towers to be exposed to the cleaning solution for a contact time of 10 minutes before powering on the NEO Iris. The solution will then automatically drain when the NEO Iris is powered on.



**Note:** After powering on the NEO Iris, the reader lamp must be allowed to warm up for at least thirty (30) minutes.

### Images

The following images illustrate this task:



Do not add solution to overflow troughs

### **Cleaning the Common Waste Container**

### Purpose

The purpose of cleaning the common waste container is to reduce gross contamination in the container.

### Required Tool(s)

The following tool(s) are required to complete this task:

• 0.25% sodium hypochlorite (NaClO) as the recommended cleaning solution

*Refer to Monthly Maintenance- Decontaminate Tubings for the specifications and preparation of 0.25% NaClO* 

• Alcohol wipes, if available

#### Procedure

Note: Only follow this procedure when the NEO Iris is switched off.

To clean the common waste container:

Step	Action
1	Prepare a concentration of 0.25% NaClO working solution.
2	Empty the common waste container.
3	Vigorously rinse the common waste container with the working solution of the recommended cleaning solution. Allow the common waste container to be exposed to the cleaning solution for a contact time of 10 minutes.
4	Use the alcohol wipes or the recommended cleaning solution to clean the inside of the common waste container lid. Allow the container lid to air dry after using the alcohol wipes, or wipe the container lid with deionized water following use of the recommended cleaning solution.
5	Empty the recommended cleaning solution out of the common waste container after the 10 minute contact time has expired.
6	Rinse the common waste container with deionized water.
7	Empty the deionized water out of the common waste container.

8	Reconnect the common waste container back onto the instrument.

## Inspecting the Syringes

### Purpose

The purpose of inspecting the syringes is to determine if any of the syringes are too loose and/or are leaking.

### Procedure



Note: Only follow this procedure when the NEO Iris is switched off.

To inspect the syringes:

Step	Action	
1	Inspect the syringes for leakage. Syringes must not be leaking.	
2	Check the barrels of the syringes to determine if they are finger tight by grasping the knurled screw and twisting gently clockwise. You must never use tools for this procedure. Syringes must be finger-tight secured to the bottom of the 3-way valve.	
Close-up of knurled securing screw at top of syringe connected to 3-way valve		
	Note: If you determine that, for whatever reason, a syringe must be replaced, you r follow the <b>Removing and Replacing a Syringe</b> procedure described later in this ch perform the necessary follow-up maintenance tasks as described in that procedure.	

# Monthly Maintenance Tasks

# Monthly Tasks

The following tasks are to be completed on a monthly basis:

- Clearing Test Data
- Decontaminate Tubings

## **Clearing Test Data**

#### Purpose

The purpose of deleting statistical test data via the Utilities>Statistics tab is to delete data that is not deleted using the archive method. If such data is allowed to build up, this could impact the PC performance. The statistical data is related to the Tests per Assay report. This clearing task should be performed monthly, but not until after your site statistical analysis period has ended. If your analysis period is greater than a month, then the frequency of this task will be dictated by that analysis period schedule. Once data is cleared, it will no longer be available for inclusion in a future Test per Assay report. Data conservation should be based on your analysis period schedule. If this data is not required to be analyzed by your site at any time, then this clearing task should be performed monthly.

### Procedure

To clear test data:

Step	Action
1	Press the <b>Utilities</b> button on the Main Menu Bar.
2	Select the <b>Statistics</b> tab.
3	Select the <b>Tests per Assay</b> report name. The <b>Clear Test Data</b> button will then be displayed.
4	Press the <b>Clear Test Data</b> button. A confirmation dialog will then be displayed.
	Clear Test Data
5	Press the <b>Yes</b> button of the confirmation dialog to proceed or the <b>No</b> button to cancel the
	task.
	Clear Test Data?         It will clear data relating to the number of tests performed per assay within the date ranges selected.         It will not affect result data, but reports of this type produced in the future will not include all values.         Are you sure you wish to delete this data?
	Yes No

Step	Action
6	Press the <b>OK</b> button of the <b>Test data cleared</b> message to acknowledge the completion of the task.
7	Record completion of the task on the monthly maintenance record. Note: If the data is not cleared on a monthly basis, based on the requirements of y statistical analysis period, then record N/A (or whatever equivalent nomenclature yo documentation requires) on the monthly maintenance record to indicate that it was performed within a given month. However, you may perform this task on a subseque month, as defined by your analysis period. N/A = Not Applicable.

### **Decontaminate Tubings**

#### Purpose

The purpose of decontaminating tubings is to meet biohazard regulations.

### Required Tool(s)

The following tool(s) are required to complete this task:

- 0.25% sodium hypochlorite (NaClO) as the recommended cleaning solution
- PBS
- 1L container
- Alcohol wipes, if available
- WipeAll or paper towels

<u>Note</u>: You must still keep your original NEO Iris 10 L system liquid container available, even If you use cubed saline as your system liquid reservoir, because that original NEO container will be required for the decontamination procedure.

### Specifications for Sodium Hypochlorite (NaCIO)

NaCIO sourced for use in the decontamination protocol should be aligned to the following specifications:

CAS Number:	Sodium hypochlorite: 7681-52-9
Grade	Commercial
Concentrations	Sodium hypochlorite: Variable



**Warning:** Common household bleach or sodium hypochlorite with detergent additives should not be used as these could affect the performance of the system.

### Preparation of 0.25% Sodium Hypochlorite (NaCIO) Solution

- 2 L of 0.25% NaClO is required for monthly maintenance
- The dilution is prepared using deionized water

- The concentration for the purchased NaCIO solution should be verified on the certificate of analysis (CoA) for the lot used
- To determine the volume of NaCIO required based on the concentration in the CoA, use the formula (V1) = (C2)(V2)/(C1) where
  - V1 = volume of undiluted NaClO required
  - C2 = concentration of NaCIO desired (0.25%)
  - V2 = volume of 0.25% NaClO desired (2000mL)
  - $\circ$  C1 = concentration of NaCIO indicated on the CoA
  - o Otherwise, V1(mL) = 500/C1

Volumes required for typical NaClO			
C1 NaClO Concentration	DiH2O	V1 Volume Undiluted NaClO	Final Volume of Prepared Solution
8.27%	1940 mL	60 mL	
12.0%	1958mL	42mL	2L
14.0%	1965mL	35mL	



**Warning:** Washer dispense check failures may occur during decontamination maintenance checks. If the ion concentration of a solution is too low, the liquid level detection sensors may fail to detect the presence of decontamination solution used on the instrument. Should you observe repeated washer dispense check failures after verifying the NaClO dilution, and have ruled out other mechanical related causes for the error, it may be necessary to perform dilution of the NaClO using phosphate-buffered saline (PBS) instead of DI water.

### Procedure

To decontaminate the tubings:

Step
------

Action

Step	Action
1	Press the <b>Maintenance</b> button on the Main Menu Bar.
2	Select <b>Decontamination tubings</b> from the action list.
3	Press the <b>Start</b> button.
4	Disconnect the blue tubing connected to the 20 L system liquid container and place it into the empty 1 L collection container. This is the return tube from the de-bubbler mechanism.
5	Put a plate into the plate loading tower, manually enter barcode identification data, and assign the assay <b>DB_CK</b> to it. Start the assay.
6	After the <b>DB_CK</b> assay has completed and the plate has been returned to the plate loading tower, verify that liquid has collected into the empty 1 L container from step 4. Record this verification on the monthly maintenance record.
7	Empty the liquid that has collected into the 1 L container from the blue tubing. Note: If there is no liquid in this 1 L container after the DB_CK assay has completed and it has been verified the blue tubing was correctly placed inside the container, then contact Technical Support. This may indicate a failure of the de-bubbler mechanism.
8	Connect the blue tubing back to the 20 L system liquid container lid.

Step	Action
9	Disconnect the system liquid containers at the two (2) in-line connectors (circled in the photograph below) and store the disconnected in-line filter in a secure place ready for reinstallation.
	System liquid containers connected
	System liquid containers disconnected
1	Disconnected in-line filter
	Note: If any part of the in-line filter or its associated tubing is discolored it should be
	discarded and a new filter and/or its associated tubing should be used when it is reinstalled.

Step	Action
10	Remove the cap and tubing from the 20 L system liquid container and set aside on a clean paper towel. Remove the 10 L and 20 L system liquid containers from the system. Empty the PBS from both containers and fill the 20 L container with 2 L of the recommended cleaning solution.
11	Swish the cleaning solution vigorously in the 20 L container, ensuring all interior areas of the container are coated. Tilt the container forward to allow the solution to fill the inline filter connection port.

Step	Action
12	Return the 10 L system liquid container to the NEO Iris. Place the system liquid cap and tubing into the 10 L container for decontamination. All tubing and the sensor post should be in the container. It is not necessary to screw the lid onto the container.  Note: When placing the cap and tubing into the 10 L container the sensor may indicate insufficient volume; however there will be enough decontamination solution in the container to allow the assay to run and avoid additional alarms. Any system liquid empty alarms should be acknowledged.  20 L container lid and tubing placed in 10 L container
13	Put a plate into the plate loading tower, manually enter barcode identification data, and assign the assay <b>PrimeDec</b> to it. Start the assay.
14	The instrument will aspirate the recommended cleaning solution and allow it to soak within the system tubing for the 10 minute recommended soak time. The <b>PrimeDec</b> assay will complete once the recommended soak time has expired, and the plate will return to the plate loading tower

Step	Action
15	The 10 L and 20 L system liquid containers can each be rinsed when the plate returns to the plate loading tower at the completion of the PrimeDec assay. Remove the cap and tubing from the 10 L system liquid container and allow the solution to drain from the tubing back into the container. Set the cap and tubing aside on a clean WipeAll or paper towel.
	Rinse both containers with a minimum of 1000 mL of PBS a minimum of two (2) times. <b>Tap</b> both 10 L and 20 L containers to ensure any residual cleaning solution is removed from the inline filter connection port each time the containers are rinsed and emptied completely.
16	Return both of the empty containers to the instrument; however do not reconnect the system liquid containers at this time. Note: So that you can optimize the use of the cleaning solution, during the decontamination period, you can empty the common waste container and then decontaminate it with the recommended cleaning solution that was used to clean the system liquid containers. The container must then be rinsed with deionized water prior to usage. Use some alcoholl wipes or the recommended cleaning solution to clean the inside of the common waste container lid, followed by deionized water. Return the clean common waste container to its proper place on the instrument. This waste container cleaning task must be performed on a weekly basis. Refer to <b>Cleaning the Common Waste Container</b> in this chapter for more details regarding the weekly cleaning of the common waste container.

Step	Action
17	With he cap and tubing still set aside, place only the de-bubbler tubing that extends from the bottom of the cap (de-bubbler output tubing) into the empty 1 L container. Fluid will collect in this container in step 18.  Note: The sensor float must continue to be lifted to prevent system liquid empty warnings.  Sensor post and float De-bubbler output tubing
18	Put a plate into the plate loading tower, manually enter barcode identification data, and assign the assay <b>ClrTube1</b> to it. Start the assay.
19	After the <b>ClrTube1</b> assay has completed, drain the remaining cleaner from the system liquid cap tubing into the collection container. Empty the 1L collection container.
20	Wipe the system liquid sensor (post and float) and the exterior of the tubing extending below the cap with a clean alcohol wipe, or recommended cleaning solution, followed by wiping with deionized water. Allow them to dry.

Step	Action
21	Fill the 10 L container with 9 L of PBS. Place the system liquid cap and tubing in the 10 L container. Place the de-bubbler tubing that extends from the bottom of the cap (de-bubbler output tubing) into the 1 L collection container.
22	Put a plate into the plate loading tower, manually enter barcode identification data, and assign the assay <b>Prime1</b> to it. Start the assay.
23	After <b>Prime1</b> has completed, remove the cap and tubing from the 10L container. Allow any remaining liquid in the tubing to drain into the collection container. Set cap and tubing aside on a clean WipeAll or paper towel. Empty the 1 L collection container and place the debubbler tubing back into it. <u>Note:</u> The sensor float must continue to be lifted to prevent system liquid empty warnings.
24	Put a plate into the plate loading tower, manually enter barcode identification data, and assign the assay <b>ClrTube2</b> to it. Start the assay.
25	After the <b>ClrTube2</b> assay has completed, drain any remaining liquid from the system liquid cap tubing into the collection container. Empty the 1 L collection container.

Step	Action
26	Add 3L more of PBS to the 10 L container. Place the system liquid cap and tubing back into the 10 L container. Place the de-bubbler tubing that extends from the bottom of the cap (de-bubbler output tubing) into the 1 L collection container.
27	Put a plate into the plate loading tower, manually enter barcode identification data, and assign the assay <b>Prime2</b> to it. Start the assay.
28	After <b>Prime2</b> has completed, empty the collection container. Place the debubbler tubing back into the empty collection container.
29	Load a clean untreated microplate into the tower, manually enter barcode identification data, and assign the assay <b>Prime3</b> to it.
30	Load a fresh centrifuged blood sample onto the NEO Iris with the barcode turned so that the NEO Iris cannot scan it. Use the <b>Recall</b> function to manually enter the barcode identification as <b>sample</b> . The typing of the word <b>sample</b> is case sensitive (all letters must be typed as lowercase). Start the assay ( <b>Prime3</b> ).

Step	Action
31	When the <b>Prime3</b> assay completes, verify that there is no visible hemolysis in the wells of the plate. Record this verification on the monthly maintenance record. If there is visible hemolysis, the following troubleshooting techniques should be implemented:
	<ul> <li>The Prime3 assay may be repeated. If hemolysis persists, the system liquid container may be checked for residual chlorine.</li> </ul>
	<ul> <li>The PBS in the 10 L container may be checked for residual chlorine. If chlorine is present, discard PBS from 10 L, refill with fresh PBS, and repeat the <b>Prime3</b> assay. Verify no hemolysis is present before continuing.</li> </ul>
32	Once it is verified that hemolysis is not present, empty the 1 L collection container and the remaining PBS from the 10 L container. Tap the 10 L container to ensure any residual liquid remaining in the inline filter connection port is removed and emptied.
33	Fill both the 20 L and 10 L system liquid containers with new PBS and return both containers to the system.
	monthly maintenance action and not as preparation for shipping the NEO Iris.
34	Connect both of the system liquid containers together using the two (2) in-line connectors (refer to step 9 for photographs). This connecting component, bridged by the single in-line filter, was previously stored in a secure place.

Step	Action
35	Return the system liquid cap and tubing to the 20 L system liquid container (including debubbler tubing). It is recommended that the reconnected tubing (and the liquid level sensing cable) be threaded though the handle of the 20 L container to prevent any potential rubbing of the tubing on the inside casing of the cabinet as the drawer is routinely pulled into and out of the cabinet.
	Make sure that the system liquid sensor float can move freely along the post, and is in the down position before securing the system liquid cap to the 20 L container.
36	Verify that the caps are properly secured on the system liquid containers. The NEO Iris is now ready for use.

# As Needed Maintenance Tasks

### As Needed Tasks

The following tasks are to be completed on an as needed basis:

- Flushing the System Liquid
- Emptying the 20-liter System Liquid Container
- Performing the Pipettor Verification Test (PipTest)
- Performing the Positions Check (PosCheck)
- Performing the Syringe Change
- Checking the Residual Volume
- Removing and Reinserting the Manifold
- Using the Washer Teach Tool
- Checking the Manifold Probes
- Adjusting the Manifold Grub Screws
- Cleaning the Liquid Level Sensor
- Cleaning the Liquid Overflow Detection Mat
- Removing and Replacing a Syringe
- Removing and Replacing a Probe
- Cleaning the Reader Mirror and Light Diffuser
- Changing the Camera Lamp
- Replacing the Y-pusher

# Flushing the System Liquid

### Purpose

You should flush the system liquid lines with deionized water if you are shutting down the instrument for an extensive period of time.

This task can also be performed to flush deionized water out of the system liquid lines and replace with PBS.



**Note:** You can also perform the Flush System Liquid maintenance action as a total system prime to remove any air bubbles that may have entered the line. You can leave the wash buffer lines connected to the container and run PBS through the lines.

It is possible for a washer dispense error to occur while the *Flush Liquid System* action is being performed. The on-screen instructions describe this risk, but explain that this is a benign event and some additional error recovery instructions for that scenario are provided.

### Procedure

To flush the system liquid:

Step	Action
1	Press the <b>Maintenance</b> button on the Main Menu Bar.
2	Select Flush Liquid System from the action list.
3	Press the <b>Start</b> button.
4	Disconnect the clear washer and system liquid tubing, and set them aside.
5	Connect the blue tubing to the input connector on top of the waste container.
6	Place a plate in the tower, manually enter a plate ID, and assign the <b>PrimeAir</b> assay to it. Start the assay.

Step	Action
7	After the <b>PrimeAir</b> assay is completed, place the terminal portions of the clear washer and system liquid tubing into a container with at least 1.5 L of your chosen flushing liquid, e.g. de-ionized water.
	<b>Note:</b> The blue tubing should still be connected to the waste container.
8	Place a plate in the tower, manually enter a plate ID, and assign the <b>PrimeDI</b> assay to it. Start the assay.
	<b>Note:</b> After the assay completes, the system liquid throughout the NEO Iris will be y chosen flushing liquid, e.g. de-ionized water.
9	After the <b>PrimeDI</b> assay is completed, reconnect the clear washer and system liquid tubing to their original connections.
10	The blue tubing should remain connected to the waste container until the system liquid in the NEO Iris has been returned to PBS.

## **Emptying the 20-Liter System Liquid Container**

#### Purpose

The purpose of emptying the 20-Liter System Liquid Container is to prepare the container prior to replacing liquid in it.

### **Before You Begin**



**Warning:** Before you can physically empty the 20-liter system liquid container, you must remove it from the washer software. If you do not adjust the Wash Buffer dialog to reflect the buffer removal and a test is in process, the software still attempts to use the wash buffer, even when none is present. If this occurs, the associated NEO Iris tubing and pumps can run dry.

#### Procedure

To empty the 20-liter System Liquid Container:

Step	Action
•	

1	Open the Wash Buffer dialog.	
2	If the wash pump is in use, select the <b>Remove</b> button. The <b>red In use</b> circle turns green.	
3	Remove the 20-liter system liquid container and empty it according to your standard laboratory procedure.	
	The green <b>Fill Level</b> check mark (✓) for the wash buffer changes to a red exclamation symbol (●).	
4	Click <b>Done</b> to exit the Wash Buffers dialog.	
	Attention: You must prime the pump after placing the 20-liter system liquid container back into the system prior to the pump's use.	
5	Prime the NEO Iris	
	Click the Maintenance button to access the Maintenance section of the software.	
	Select Flush System Liquid.	
	Refer to <b>Flushing the System Liquid</b> for detailed instructions on how to start this maintenance action.	

## Performing the Pipettor Verification Test (PipTest)

#### Purpose

The purpose of performing the pipettor verification test is to verify the precision pipetting of each probe. This verification is required during NEO Iris installation and after a syringe exchange.



**Note:** If the PipTest fails, the operator must inspect the probe(s), tubing and syringe(s) for signs of incorrect installation and leaking connections, and correct any observed errors. The PipTest must then be repeated.

### Required Tool(s)

The following tool(s) are required to complete this task:

- Scale
- Unused test tube

### Procedure

To perform the pipettor verification test:

Step	Action	
1	Press the <b>Maintenance</b> button on the Main Menu Bar.	
2	Select <b>Pipettor Verification Test</b> from the action list.	
3	Press the <b>Start</b> button.	
4	Follow the on-screen instructions:          Pipettor Verification Test <piptest>:         1. Weigh each strip of a microplate and record each value in the PipTest chart.</piptest>	
	<ol> <li>Load the plate, with all the strips, in to the tower.</li> <li>Manually enter an id and assign the 'PipTest' assay.</li> <li>Fill a tube with 7ml DIH20, place in a sample rack and load on the instrument.</li> <li>Using the Recall Function, manually enter 'tube1' as the id.</li> <li>Start the assay.</li> <li>When the assay completes, remove the microplate and weigh each strip.</li> <li>Record the weight of each strip in the PipTest chart.</li> <li>The difference of these values is the total weight. This value should be within the expected range for the strip number.</li> </ol>	

*	Extra instructions for PipTest Sample Loading (steps 4-5)	
	Step	Instructions
	4	Load a tube with 7 ml of de-ionized water on the instrument with the barcode turned so that the system cannot scan it.
	5	By using the Recall function, manually type <b>tube1</b> as the ID in the correct position of the rack. The typing of the word <b>tube1</b> is case sensitive (all letters must be typed as lowercase).

### **Recording PipTest Results According to Dispensing Pattern**

The PipTest chart is designed to record results according to the instrument dispensing pattern, as illustrated below.



# Performing the Positions Check (PosCheck)

### Purpose

The purpose of performing the positions check is to verify the probes can access all the required areas within the 14-lane Bay, 5-lane Bay, plate pipetting positions, and wash stations. This verification is required during NEO Iris installation and after a probe exchange.

### Procedure

To perform the pipettor positions check:

Step	Action	
1	Press the <b>Maintenance</b> button on the Main Menu Bar.	
2	Select <b>Positions Check</b> from the action list.	
3	Press the <b>Start</b> button.	
4	Follow the on-screen instructions:	
	<ol> <li>Load a sample rack in lane 1 and name the first position L1 and the last position L16 using the recall function.</li> <li>Load a sample rack in lane 14 and name the first position R1 and the last R16 using the recall function.</li> <li>Load the 5-position reagent rack in the far left lanes of the reagent bay.</li> <li>Name the first position LR1 and the last position LR5 using the recall function.</li> <li>After loading the rack back on the instrument enter a volume of 1ml for these 2 positions.</li> <li>Load the 9-position reagent rack in the far right lanes of the reagent bay.</li> <li>Name the first position RR1 and the last position RR9 using the recall function.</li> <li>After loading the rack back on the instrument enter a volume of 1ml for these 2 positions.</li> <li>Load the 9-position RR1 and the last position RR9 using the recall function.</li> <li>Name the first position RR1 and the last position RR9 using the recall function.</li> <li>After loading the rack back on the instrument enter a volume of 1ml for these 2 positions.</li> <li>Place a plate in the tower and manually identify it then assign the test PosCheck to it.</li> <li>Confirm the assay finishes without error.</li> </ol>	

# Performing the Syringe Change

### Purpose

The purpose of performing the syringe change is to have the instrument move the syringe drive down on all diluter pumps so a syringe can be replaced if necessary.

### Procedure

To perform the syringe change:

Step	Action
1	Press the <b>Maintenance</b> button on the Main Menu Bar.
2	Select Syringe Change from the action list.
3	Press the <b>Start</b> button.
4	Follow the on-screen instructions:
	Syringe Change <pumps down="">:         1. Load a microplate in the tower and manually enter an id.         2. Assign the 'SyringeEx' assay and start the test.         3. After the assay is finished, turn off the instrument.         4. The syringe can now be exchanged using the procedure found in the Maintenance chapter of the Operator's Manual.</pumps>
5	Follow the procedure for <b>Removing and Replacing a Syringe</b> that is included in this chapter.

### Checking the Residual Volume

#### Purpose

Checking the residual volume is an as needed maintenance action used for troubleshooting purposes, such as when unexpected results are obtained. The purpose of checking the residual volume is to verify that the volume of liquid left in the wells is correct after the wash sequence.



Note: You must only perform this maintenance action at the request of Technical Support.

### Required Tool(s)

The following tool(s) are required to complete this task:

- Scale
- Capture-R® Ready-Screen® or Capture-R<sup>®</sup> Select plate (see notes below)



**Note:** Either expired or in-date Capture-R<sup>®</sup> Ready-Screen<sup>®</sup> plates can be used. The same plate can be re-used up to ten (10) times if the plate dries between usages.

<u>Using a wet plate for multiple runs during a single troubleshooting session</u>: The same plate can be re-used up to twenty (20) times if the plate remains wet between usages (this could occur if there are failures that require adjustments and residual volume is run numerous times in a single day). The plate can either be tapped out between runs or re-run with the previous residual volume in the wells. After completion of a single troubleshooting session using a wet plate multiple times, the wet plate must be discarded.



**Note:** In addition to Capture-R<sup>®</sup> Ready-Screen<sup>®</sup> plates, you are permitted to use Capture-R<sup>®</sup> Select plates to perform *Residual Volume*. However, unlike the Capture-R<sup>®</sup> Ready-Screen<sup>®</sup> plates, if a Capture-R<sup>®</sup> Select plate is used once, it must be discarded after that single use and a new plate must be used if the maintenance action needs to be repeated. Multiple use of the Capture-R<sup>®</sup> Select plates is not permitted.

#### Plate types other than Capture-R<sup>®</sup> Ready-Screen<sup>®</sup> or Capture-R<sup>®</sup> Select are not acceptable for use

#### Procedure

To check the Residual Volume:

Step	Action
------	--------
Step	Action
------	---
1	Press the <b>Maintenance</b> button on the Main Menu Bar.
2	Select <b>Residual Volume</b> from the action list.
3	Press the <b>Start</b> button. The on-screen instructions will be displayed.
4	Press the <b>Continue</b> button. Select <b>ResVol</b> in the <i>Resource Overview</i> window. The system indicates that a plate is required to run this assay.
	Resource Dverview       Plates       Reagents       Controls       Donors       Washbuffer       Pipettor       Incubator         Image: ResVol       full       <
5	Press the <b>Plates</b> button above the exclamation mark. The system displays the <i>Plate Loading</i> <i>Tower</i> dialog.
6	Weigh a Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> plate with the scale and then correctly place the plate into a transport frame. Record the weight on the <b>NEO Iris Residual Volume Calculation Record</b> .
	Note: If a wet plate is being re-used for multiple runs during a single troubleshooti session, you should use the original dry plate weight on the calculation record (line each of the wet run calculations during the session.
7	Load the plate in any position of the Plate Loading Tower.
8	Press the Assay Selection tab and then press the Plate ID field next to the orange LED.

Step	Action
9	Type a unique plate ID, such as current date, and choose the <b>ResVol</b> assay button from the assay selection area to the right. The system displays the selected assay next to the plate ID under the Assay Code column.
	Scan meripides     No Status     Plate ID     Assay Code       15     14       13     RESVOL050009       13     RESVOL050009       11     10       10     10       9     11       10     10       7     10       6     11       2     11       10     11       10     11       10     11       10     11       10     11       10     11       11     11       12     11       10     11       11     11       12     11       10     11       11     11       12     11       10     11       11     11       12     11       13     11       14     11       15     11       16     11       17     11       18     11       19     11       10     11       10     11       11     11       12     11       13     11       14     11       15     11
10	Press the <b>Done</b> button. The system redisplays the <i>Resource Overview</i> window indicating that resources are sufficient to run the assay.
	Resource Dverview       Resources     Assay Name     Samples/Strips     Plates     Reagents     Controls     Donors     Washbuffer     Pipettor     Incubator       Res/Vol     full     Image: Controls     Donors     Image: Controls     Image: Controls<
11	Select <b>Start</b> . The transport system takes the plate to the washer module. The plate is returned to the tower when the task is finished.
12	Remove the plate from the tower and weigh it on the scale. Record the weight on the <b>NEO</b> Iris Residual Volume Calculation Record.
13	Subtract the weight from step <b>7</b> from the weight from step <b>13</b> . Record this calculated weight on the residual volume calculation record. Report this calculated weight to Technical Support.

Step	Action
14	Remove strips 1, 6, and 12 from the plate to verify that the remaining volume is even by holding the strips parallel to each other.
	If the volume is uneven, verify that the manifold is correctly mounted. If the problem persists, contact Technical Support to determine further action.
	<b>Note:</b> The residual volume is out of range if the residual volume between any two wells (within a strip and between strips) doubles (or greater).
	Excess liquid in test
	Note: If the residual volume verification fails, Technical Support may recommend that the manifold must be removed so that the manifold probes can be checked for correct positioning using the manifold align tool. If this course of action is recommended by Technical Support, all of the prescribed maintenance tasks must be subsequently performed once the manifold is reinserted back into NEO Iris, which includes repeating the residual volume verification.
	Warning: The main maintenance screen indicates that this action is Due if the assay successfully completes. The system does not analyze the results on the instrument, so you must determine if the volume is acceptable. You are responsible for making the corrections necessary to satisfy these requirements.
L	1

## **Removing and Reinserting the Manifold**

#### Purpose

The manifold can be removed and reinserted so that it can be flushed when a blocked manifold needle is suspected or when the Flow350 and/or the Flow200 maintenance tasks fail.

#### Phases

This procedure consists of two phases:

- Removing and Cleaning the Manifold
- Reinserting the Manifold

#### **Required Materials**

The following materials are required to complete this task.

- Absorbent paper towel
- Stylus
- Plastic syringe
- Probe check and align tool

#### Removing and Cleaning the Manifold

To remove and clean the manifold:

Step	Action
1	Turn off the NEO Iris module.
2	Place an absorbent paper towel in the liquid overflow detection mat tray to absorb any drops of fluid that may fall when the manifold is moved to a position above the tray in the following step 3.

Step	Action
Step 3	Action If the manifold probes are located down inside the trough, lift the manifold up above the trough to give the probes clearance and then gently slide the manifold to the right along the two (2) parallel x-rails so that is it removed from the direct area of the trough. If the trough is a start of the trough is a
	Note: Attempting to remove the manifold with the probes located down inside the trough without first lifting the probes clear and away from the trough can result in bent probes
	tough can result in bent probes.

Step	Action
4	Manually loosen the thumb screw used to secure the manifold to its support arm by turning it counter-clockwise with your fingers and thumb until it detaches from the support arm. The location of the thumb screw is shown in the photograph below.
	<b>Note:</b> The thumb screw will stay attached to the manifold front handle by the use of a containment device on the reverse side of the handle.
	Note: Do not adjust the two (2) securing screws located at the bottom of the manifold handle in any way. The location of the two (2) securing screws is shown in the photograph below.
	0
5	Grip the top handle tab and pull the manifold towards you. The manifold will detach from a small metal post. A small gray post will also be visible at the back of the support arm.
	Metal post and small gray post visible on manifold supporting arm
6	Remove and discard the absorbent paper towel that was placed in the liquid overflow
1	detection mat tray to absorb any drops of fluid in step 2.

Step	Action
7	To flush the manifold with PBS to unblock potentially blocked probes, you must use a plastic syringe filled with PBS to push PBS through the holes in the top of the manifold that supply the dispense probes and carry away the fluid from the aspiration probes. A stylus can be used to stylus out any resistant obstructions in the probes that cannot be removed by flushing with PBS. Two (2) sizes of stylus are provided, one for each probe type.
	9
	Stylus
	Stylus the manifold
	Flushing the manifold

Step	Action
8	Use the probe check and align tool to check the manifold probes. Straighten probe(s) if needed by gently pushing the necessary probe to a vertical position.

# Reinserting the Manifold

To reinsert the manifold:

Step	Action
1	Grip the top handle tab and push the manifold away from you into the manifold support arm so that the guiding tips on either side at the back of the manifold correctly engage the holding assembly.
	Note:Make sure that the manifold engages the small metal post and pushesback against the small gray post at the rear of the assembly.
2	Manually tighten the thumb screw to secure the manifold to its support arm by turning it clockwise with your fingers and thumb until it tightens into the support arm. The location of the thumb screw is shown in the photograph below.
3	Power on the NEO Iris module.
4	Use the washer teach tool to teach the washer manifold its zero position.
5	Perform the weekly Verifying the Washer maintenance task.
6	Perform the <b>Checking the Residual Volume</b> maintenance task.

## Using the Washer Teach Tool

#### Purpose

The purpose of using the washer teach tool is to teach the washer manifold its zero position. This task must be performed on an as needed basis, when the manifold has been removed and replaced or when the **Checking the Residual Volume** maintenance task fails.

### Required Tool(s)

The following tool(s) are required to complete this task:

• 2 mm thick washer teach tool



#### Procedure

To use the washer teach tool:

Step	Action
1	Press the Maintenance button on the Main Menu Bar.
2	Select <b>Teach Washer</b> from the action list.
3	Press the <b>Start</b> button. The onscreen instructions are displayed.
4	Press the <b>Continue</b> button. The <b>Resource Overview</b> window is displayed.

Step	Action
5	Select the <b>TeachWash</b> assay, which will highlight in <b>blue</b> , and press the <b>Plates</b> button above the <b>red</b> exclamation mark. The <b>Plate Loading Tower</b> dialog will be displayed.
6	Place the 2 mm thick washer teach tool fully down into a plate carrier, with the <b>red</b> dots facing DOWN.
7	Place the loaded plate carrier into the plate loading tower into any available tower position.
8	Press the <b>Assay Selection</b> tab and select the <b>Plate ID</b> field next to the <b>orange</b> LED onscreen signal.
	Plate Loading Tower     No Status     Plate ID     Assay Code     Strip Selection     Scan Plate     Assay Selection     Re       Schedule     13     Image: Schedule     Flow 200     Flow 350     Re
9	Type a unique plate ID, such as the current date, and assign the <b>TeachWash</b> assay from the buttons to the right. The system displays the selected assay next to the <b>Plate ID</b> under <b>Assay Code</b> .
10	Press the <b>Done</b> button on the <b>Plate Loading Tower</b> dialog. The system returns to the <b>Resource Overview</b> window indicating that all resources are sufficient to run the assay.

Step	Action
11	Press the <b>Start</b> button on the <b>Resource Overview</b> window. The transport system takes the teach tool to the washer module. The manifold moves down to the teach tool at programmed positions. When the manifold detects the tool, that position is the value the instrument teaches as zero. This value is used when the manifold aspirates liquid form the wells. At the end of this process, the tool is returned back to the plate tower.           Mote:         After the successful completion of the TeachWash assay, the main maintenance status screen indicates that Teach Washer is <i>Due</i> . The state remains unchanged.
12	Remove the tool from the plate tower and then the transport frame, and store it in a secure place for future use.

## **Checking the Manifold Probes**

#### Purpose

The purpose of checking the manifold probes is to ensure that the probes are correctly positioned and not bent. This task must be performed on an as needed basis, when the manifold has been removed and replaced or when the **Checking the Residual Volume** maintenance task fails.

### Required Tool(s)

The following tool(s) are required to complete this task:

• Probe check and align tool





## Procedure

To use the probe check and align tool:

Step	Action
1	Once the washer manifold is removed from the NEO Iris, slide it down onto the probe check and align tool.
	Note: The manifold will only securely slide down onto the tool in one orientation. The longer set of probes will slide down into the deeper tool groove and the shorter set of probes will slide down into the shallower tool groove on the opposite side of the tool. The two (2) grub screws (one on either end of the manifold) sit flush down onto the top of the tool.

<ul> <li>Check that the manifold probes are straight. Make sure that:         <ul> <li>All probe tips should be flush with the respective surfaces of the tool.</li> <li>The probes should be parallel to one another.</li> </ul> </li> <li>The tool defines the spacing between aspirate and dispense probes.</li> <li>Slight corrections, by pushing the probes, are permitted to straighten any bent or misaling probes.</li> <li>If probes are kinked, then the manifold must be replaced.</li> </ul>	Step	Action
<ul> <li>a. All probe tips should be flush with the respective surfaces of the tool.</li> <li>b. The probes should be parallel to one another.</li> <li>The tool defines the spacing between aspirate and dispense probes.</li> <li>Slight corrections, by pushing the probes, are permitted to straighten any bent or misaliprobes.</li> <li>If probes are kinked, then the manifold must be replaced.</li> </ul>	2	Check that the manifold probes are straight. Make sure that:
<ul> <li>b. The probes should be parallel to one another.</li> <li>The tool defines the spacing between aspirate and dispense probes.</li> <li>Slight corrections, by pushing the probes, are permitted to straighten any bent or misaling probes.</li> <li>If probes are kinked, then the manifold must be replaced.</li> </ul>		a. All probe tips should be flush with the respective surfaces of the tool.
The tool defines the spacing between aspirate and dispense probes. Slight corrections, by pushing the probes, are permitted to straighten any bent or misaling probes. If probes are kinked, then the manifold must be replaced.		b. The probes should be parallel to one another.
Slight corrections, by pushing the probes, are permitted to straighten any bent or misaling probes. If probes are kinked, then the manifold must be replaced.		The tool defines the spacing between aspirate and dispense probes.
If probes are kinked, then the manifold must be replaced.		Slight corrections, by pushing the probes, are permitted to straighten any bent or misaligned probes.
		If probes are kinked, then the manifold must be replaced.
3 Remove the manifold from the tool so that the manifold can be reinserted back into the Iris washer module.	3	Remove the manifold from the tool so that the manifold can be reinserted back into the NEO Iris washer module.

# Adjusting the Manifold Grub Screws

#### Purpose

The purpose of adjusting the manifold grub screws is to ensure that the probes are correctly positioned in relation to the wash aspiration position when washing a plate. This task must be performed on an as needed basis when the manifold probes are checked and the probe tips are found not to be flush with the respective surfaces of the tool.



**Note:** This task should only be performed after express direction is received from Technical Support and under no circumstances should you undertake this task without that direction.

### **Required Tool(s)**

The following tool(s) are required to complete this task:

- Probe check and align tool
- 1.5 mm Allen wrench



## Procedure

To adjust the manifold grub screws:

Step	Action
1	Once the washer manifold is removed from the NEO Iris, slide it down onto the probe check and align tool.
	Note: The manifold will only securely slide down onto the tool in one orientation. The longer set of probes will slide down into the deeper tool groove and the shorter set of probes will slide down into the shallower tool groove on the opposite side of the tool. The two (2) grub screws (one on either end of the manifold) sit flush down onto the top of the tool.

Step	Action
2	If probe tips are not flush with the respective surfaces of the tool, use the 1.5 mm Allen wrench to adjust one or both of the grub screws as needed. The wrench slots in through the holes at each end of the tool and then engages with the grub screw directly adjacent to the hole. Once engaged, the wrench can be used to adjust the grub screw.
	A longer grub screw leads to more residual volume because the aspiration probes are raised relative to the plate surface.
	Note: Only make adjustments of less than a quarter (¼) turn at any one time.
3	Remove the manifold from the tool so that the manifold can be reinserted back into the NEO Iris washer module.
4	Once the NEO Iris is ready for use, run the residual volume maintenance task. If the residual volume results are unacceptable, remove the washer manifold once again and repeat steps 1, 2 and 3 of this procedure until an acceptable residual volume result is obtained.

# **Cleaning the Liquid Level Sensor and Trough**

#### Purpose

The Liquid Level Sensor (LLS) can be removed from the washer manifold trough, cleaned and reinserted if errors are generated that indicate that the washer manifold has aspiration or dispensing check failures, or that the LLS board is reporting an error. The LLS is used to measure conductivity in the washer manifold trough so that excess of liquid in the trough can be detected in the event of manifold aspiration or dispensing problems. There are eight (8) pairs of short and long fingers. Conductivity is measured between the pad at the end of the long finger (the emitter) and the two pads at the front and rear of the short finger (the two collectors). Both collector pads are at different relative levels (upper and lower) on the opposite sides of the short fingers. The LLS is one piece.





Note: This task must only be performed as needed.

### Required Tool(s) and Materials

The following tools and materials are required to complete this task.

- Allen wrench or hex key for 2 mm hex screws.
- Source of room temperature de-ionized water.
- Disposable (soft bristle only) non-electric manual toothbrush.
- Disposable cotton swab



**Note:** Tap water must not be used; to reduce the risk of the use of chlorinated water which may also contain other unidentified contaminants.

#### Phases

This procedure consists of two phases:

- Removing and cleaning the LLS
- Reinserting the LLS

### Removing and Cleaning the LLS

To remove and clean the LLS:

Step	Action
1	Turn off the NEO Iris.
2	If the manifold probes are located down inside the trough, lift the manifold up above the trough to give the probes clearance and then gently slide the manifold to the right along the two (2) parallel x-rails so that is it removed from the direct area of the trough.
	Note: Attempting to slide the manifold to the right along the two (2) parallel x-rails with the manifold probes still located down inside the trough, without first lifting the manifold probes clear, can result in bent probes.

Step	Action
3	Pull the trough towards you by using your finger to horizontally pull the trough forward from the back so that the anchoring screws uncouple from the anchor slots. Raise the trough up and away from the anchoring base.

Step	Action
4	On a flat surface away from the instrument, unscrew the two (2) screws of the LLS using the Allen wrench or hex key for 2 mm hex screws. Place the two (2) screws in a safe place ready for future use.
	Note: When handling the manifold trough, do not adjust the two (2) anchoring pin screws on the base of the trough. These two (2) anchoring pin screws are positioned to optimally allow the removal and subsequent insertion of the trough and also to ensure maximum stability when the trough is in place in the washer module.
	Front     Back       anchorin     pin screw
	screw

Step	Action
5	Clean the LLS fingers with room temperature de-ionized water using the disposable toothbrush.
	Note: Do not brush, or bring the water into contact with, the connector (located on both sides of the LLS) or the electrical board at the end of the LLS. This can cause irreparable damage to the LLS which will then require replacement.
	Connector Electrical
	The connector (located on both sides of the LLS) and electrical board are circled on the photographs below.
6	Inspect the golden sensor pads for corrosion or dirt. Repeat the LLS cleaning if necessary. Contact Technical Support if the LLS is corroded or cannot be cleaned, and therefore subsequently requires replacement.
7	Use the disposable cotton swab to clean the trough.

## Reinserting the LLS

To reinsert the LLS:

Step	Action
1	Use the 2 mm Allen wrench or hex key to loosely attach the LLS back into the trough by only partially screwing down the two (2) LLS screws, so that the screws do not fall off the LLS/trough combination unit.
	Note: This is a preparatory step so that the long LLS sensor fingers can be correctly positioned before the two (2) LLS screws are tightened down securely using the 2 mm Allen wrench or hex key in step 2 below.
2	Push the LLS to the back of the trough (right side in picture) so that the long sensor fingers touch the well side walls. The purpose is to set the shorter fingers as far away from the left side walls as possible. Tighten the two (2) screws using the Allen wrench or hex key for 2 mm hex screws while pushing the board to the side.

Step	Action
3	Rinse the priming trough wells and LLS attached to it (as shown below) with a solution of 70% isopropanol.
	Remove the residual isopropanol and rinse with DI water. Remove the residual liquid.
	Note: Only rinse parts of the sensors that are seated into the trough wells.

Step	Action
4	Hold up the end of the trough nearest you at a slight upward vertical angle so that the <b>back</b> anchoring screw and the short clear tubing are lined up with the <b>back</b> anchoring slot and the tube hole respectively in the supporting base.
	Front Back Short
	g screw screw tube
	Preparing to line up the troughBack of the trough lined up
	Lower the trough down so that both of the anchoring screws are inserted into the two (2) anchoring screw slots of the supporting base. After the trough is lined up and level, push the trough away from you, to the back, until it is locked in position (refer to the photograph below).

The series

Step	Action
5	Power on the instrument and allow Initialization to take place.
	<b>Note:</b> Repeat the LLS cleaning procedure if the Initialization generates a washer manifold error. Contact Technical Support if the error persists after this second cleaning and Initialization.

## **Cleaning the Liquid Overflow Detection Mat**

#### Purpose

The liquid overflow detection mat is orange in color and is located in the washer module. It is used to detect washer manifold fluid leaks over the plate washing area. If such leaks do occur, then a washer liquid overflow error message is displayed and the instrument is halted.



## **Required Materials**

The following materials are required to complete this task.

• Absorbent paper towel

#### Procedure

To clean the liquid overflow detection mat:

Step	Action
1	After the washer liquid overflow error message is displayed, confirm that the screen indicates that the NEO Iris is stopped by the display of <b>Halted</b> in the status bar.
2	Reach into the washer overflow mat (via the left side of the plate loading tower and above the centrifuge loading unit), with an absorbent paper towel to dry the mat completely so that there is no longer any liquid present.
3	Press the <b>OK</b> button of the error message to enter error recovery. Troubleshoot the fluid leak. Refer to <b>Chapter 11 – Troubleshooting the NEO Iris</b> for troubleshooting details.

# Removing and Replacing a Syringe

### Purpose

A syringe should be replaced when the pipettor verification test consistently fails.

## Procedure

To remove and replace a syringe:

Step	Action
1	Run the <b>Syringe Change</b> maintenance action to prepare the position of the syringes so they are ready for removal.
2	After the assay is finished, turn off the NEO Iris power.
3	Disconnect the power cord from the NEO Iris next to the power switch on the side of the instrument.
4	Open the hood and move the right pipettor arm all the way to the left side of the loading bay taking care that the probes do not hit the sides of the bay. Open access will be available to the syringes.

Step	Action
5	Place some absorbent material, such as a paper towel, under the syringe to be removed to catch a small amount of system fluid. The diagram below indicates which syringe is associated with which probe.
	Sample Probes Tagent Trobe
	Probe 1
6	If the probe of the syringe being replaced is not over a wash tower, place some absorbent material under it when the plunger is being pushed into the syringe, to collect any expressed fluid.

Step	Action
7	Pull the hinged square steel clip (on the clip assembly) of the syringe in question from the center to a left position so that the tension is released (refer to the sequence of photographs below). No tool is required for this clip release.
	Example of syringe fully locked into dilutor pump with plunger extended
	Close-up of syringe fully locked into dilutor pump using the clip assembly
	Close-up of hinged square steel clip with tension released

Step	Action
8	Open the hinged black bar from the flat position to the perpendicular position, so that it points towards you.
9	Grasp the knurled securing screw and unscrew the syringe barrel in a counter-clockwise direction to disconnect the barrel from the 3-way valve.
10	Remove the syringe assembly and set it aside.
11	Place the new syringe in position at the 3-way valve, grasp the knurled screw, and twist the syringe clockwise to secure it to the valve. You must never use tools for this procedure.
	Warning:It is critical the syringe is held in perfect vertical orientation when twisting into the 3-way valve. Do not force the syringe. If there is resistance back it off and try again.
	Warning: Verify that you are replacing the syringe with another syringe of the same size. Replacing the syringe with a wrong-sized syringe causes pipettor inaccuracy and incorrect test results.

Step	Action
12	Grasp the plunger and pull it down low enough so that the groove on the terminal metal bulb rests over the matching ridge of the dilutor pump clip assembly.
13	Close the hinged black bar back down to the flat position. The black bar should be flush with its counterpart on the diluter pump assembly (see step 12 above).
	<b>Note:</b> If the tip of the plunger is not seated correctly then the hinged black bar will not close to its flush position.
14	Close the hinged square steel clip to the right so that the hinged black bar over the syringe plunger is secured. The tip of the terminal metal bulb on the plunger must be visible below the bottom level of the hinged black bar.
	Note: If the tip of the terminal metal bulb on the plunger is not visible below the layel of the binged black bar, this can indicate that the plunger is not
	seated correctly in the clip assembly.
15	Close the hood, plug the instrument power cord back in, and turn on the instrument. Allow the instrument to fully initialize.

Step	Action
16	Perform the Pipettor Verification Test (PipTest) to confirm proper syringe performance.

## **Removing and Replacing a Probe**

#### Purpose

A probe should only be replaced when directed by Technical Support. Such scenarios can include:

- When a bent probe cannot be straightened by gentle manual bending.
- As part of an investigation of unexpected sample reactions.
- When probe adaptor problems present themselves, as diagnosed by Technical Support.

#### Required Tool(s)

The following tools are required to complete this task:

- 2 mm Allen wrench
- Pair of laboratory scissors

#### Procedure

To remove and replace a probe:

Step	Action
1	Switch off the NEO Iris module.
2	Lift up the instrument hood.
3	Move the right pipettor arm all the way to the right side of the NEO Iris module, taking care that the single probe does not hit the sides of the loading bay. Open access will be available to the syringes.


Step	Action
7	Reach into the body of the NEO Iris module and use the pair of laboratory scissors to clean- cut the flared end of the tubing off the length of the tubing.
	Note: The pair of laboratory scissors is required to cut the flared end of the tubing off because the next steps involve the pulling of this narrow tubing through extremely tight fitting sleeve tubing. If the flared end is not cut off using a pair of laboratory scissors, then the narrow tubing will not easily be pulled through the sleeve tubing as it should be able to be pulled.
8	Move the left arm to the right so that it is over the 14-lane bay, so that all the required probe adaptors can easily be accessed.
	If the probe is to be changed on the right arm, move the right arm so that the single probe is located over the 14-lane bay.
9	Spread the probes on the left arm as necessary to easily access the required probe.
10	Gently detach the white ribbon cable from the slot on the required probe assembly by pulling downwards on the white cable. The other end of the cable will remain attached to the probe adaptor. The detached end of the cable has a small yellow tab affixed.
	White ribbon cable attached toWhite ribbon cable detached fromprobe assembly slotthe probe assembly slot
11	Gently pull the probe downwards so that the probe adaptor screw is accessible with the 2 mm Allen wrench without the interference of the white ribbon cable connector area on the probe assembly.

Step	Action	
12	Unscrew the probe adaptor from the z-rod using the 2 mm Allen wrench so that the adaptor is free from the z-rod.	
	2 mm Allen wrench engaged with the screw in the probe adaptor       Use the 2 mm Allen wrench to unscrew the probe adaptor	
13	Pull the combined probe/adaptor/tubing assembly down and away from the z-rod so that t detached tubing disappears from the syringe area and the tubing is gradually pulled throug the protective sleeve tubing. Keep pulling until the probe/adaptor/tubing assembly is completely detached from the NEO Iris module. The tubing can be pulled from either direct under the z-rod or directly above the z-rod. Pulling the tubing above the z-rod may offer le resistance when compared to pulling under the z-rod.	
	<b>Note:</b> Approximately three (3) to four (4) feet of tubing will be attached to the probe adaptor when it is completely removed from the NEO Iris module.	
14	Inserting the new tubing:	
	Take the new probe and thread the end of the tubing up through the empty z-rod.	
	After the tubing is through the z-rod, push the tubing into the sleeve tubing and continue to push all the way above the z-rod. Hold the sleeve tubing located above the z-rod so that the probe tubing glides up into sleeve tubing without pushing the sleeve tubing off the z-rod.	
15	Keep pushing the tubing until the end of the tubing appears in the syringe area and the probe adaptor is flush up to the end of the z-rod.	
	Note: After the adaptor is in place on the z-rod, the guide tubing can be reseated into the z-rod if it has been dislodged from the top of the z-rod during the threading process in steps 14 and 15.	

Step	Action
16	Screw the probe adaptor onto the z-rod using the 2 mm Allen wrench so that the adaptor is secured onto the z-rod.
	2 mm Allen wrench engagedUse the 2 mm Allen wrench towith screw in probe adaptorscrew the probe adaptor onto the assembly
17	Gently attach the white ribbon cable to the probe assembly by pushing the cable upwards into the corresponding slot on the probe assembly. Hold the cable just below the yellow tab when inserting.
	<b>Note:</b> Inserting the cable with excessive force can cause kinking of the ribbon cable with possible breakage of some elements of the internal wiring.
18	Move the arm that received the new probe to its respective side of the NEO Iris module to provide unobstructed access to the syringes. Take care that the probes do not hit the sides of the loading bay.





Step	Action	
24	Switch on the NEO Iris module. Allow the instrument to fully initialize.	
	<b>Note:</b> The NEO Iris module must be switched on for at least thirty (30) minutes prior to using the camera module to read plates. The camera lamp must be provided sufficient time to warm up fully.	
25	Perform the <b>Check Pipettor Reference</b> , <b>Positions Check</b> (PosCheck) <b>Pipettor Self Check</b> and <b>Pipettor Verification Test</b> (PipTest) maintenance tasks to confirm proper probe performance.	
	Note: Look at the syringe(s) used in the probe replacement(s) to ensure that there is no leakage of fluid at the knurled screw interface on the 3-way valve. Leakage of fluid may indicate an inadequately positioned plastic washer and/or tip.	

### **Cleaning the Reader Mirror and Light Diffuser**

#### Purpose

The purpose of cleaning the reader mirror and light diffuser is to remove spills or dirt hampering the reader functionality.

You must perform this task as needed, such as when spills or dirt are identified as hampering the reader functionality.

Cleaning is also required if dirt or spills cause a test plate to be aborted at the clean read stage.

#### Required Tool(s)

The following tool(s) are required to complete this task:

- Isopropanol wipes
- Clean dry absorbent laboratory disposable wipes
- Flash light

#### Procedure



Note: Only follow this procedure when the NEO Iris is switched off.

To clean the reader mirror and light diffuser:

Step	Action
1	Power off the NEO Iris module.
2	Lift up the NEO Iris hood.
3	Move the 4-probe arm to the far right of the supporting X-rails.

Step	Action
4	<text><image/><image/><image/><table-row><table-row><table-row><table-row><table-row><table-row></table-row></table-row></table-row></table-row></table-row></table-row></text>
	top right corner of the cover plate
5	Lift the cover plate away from the body of the instrument and remove it from the instrument cavity. Store it in a safe location away from the instrument. The reader is visible when the cover is removed.

Step		Action	
6	So that the light diffuser can be expose the bottom surface of the	accessed for cleaning, the reader door	r must be opened to
	A silver rotating wheel can be us wheel is in the top position when	ed to open or close the reader door. In the door is closed (circled in yellow	The hole in the silver in the photo below).
	If the door is closed, open the data approximately 180° in a clockwis position (circled in green in the p	oor by using your finger to turn the si e direction so that the hole in the whe bhoto below).	lver wheel eel is in the bottom
	The wheel is located to the left o	f, and below, the reader door.	
	Closed position	Open position	
	Closed position Door closed direct below the lig diffuse	, y ht er	r n, w light iser
	Closed position	Open position	

Step	Action
7	With the aid of the flash light for illumination, use an isopropanol wipe to gently clean the bottom surface of the light diffuser and also the reader mirror. The reader mirror is on the bottom of the reader cavity and the bottom surface of the light diffuser is in the top of the reader cavity.         Image: Bottom surface of the light diffuser is in the top of the light diffuser is in the top of the light diffuser is in the top of the reader cavity.         Image: Bottom surface of the light diffuser is in the top of the light diffuser         Image: Bottom surface of the light diffuser         Image: Reader mirror         Reader mirror         The top surface of the light diffuser, in the bulb cavity, can also be cleaned if needed.         Image: Do not remove the light diffuser to clean it.         Image: Note: You should clean the light diffuser first, followed by the reader mirror.
8	With the aid of the flash light for illumination, use a clean dry absorbent laboratory disposable wipe to dry the bottom surface of the light diffuser and also the reader mirror (there must be no streaks left on the mirror). The top surface of the light diffuser, in the bulb cavity, can also be dried if needed.
9	Lift the cover plate onto the body of the instrument through the instrument cavity. The internal cavity of the reader is now not visible when the protective cover is in place.
10	Screw the two thumb screws clockwise to secure the protective plate in place.
11	Close the NEO Iris hood.
12	Power on the NEO Iris module.

Step	Action
13	Complete the Take Flatfield Images maintenance task.
	Limitation: The instrument must be switched on for at least thirty (30) minutes prior to beginning Take Flatfield Images maintenance task so that the reader lamp can warm up.

#### Changing the Camera Lamp

#### Purpose

The purpose of changing the camera lamp is to replace a malfunctioning lamp.

You must perform this task as needed, such as when the reader functionality is inadequate due to the lamp malfunctioning.

#### Procedure



Note: Only follow this procedure when the NEO Iris is switched off.

To change the camera lamp:

Step	Action
1	Power off the NEO Iris.
2	Lift up the NEO Iris hood.
3	Move the 4-probe arm to the far right of the supporting X-rails in the main body of the NEO Iris module.

Step	Action
4	Unscrew the two thumb screws holding the reader protective cover plate in place. Both screws are located on the top right and left corner of the cover plate respectively. These two screws are kept in place, by using internal springs to attach the screws to the cover plate when unscrewed from the body of the instrument.
	Reader protective cover plateClose-up of thumb screw in the top right corner of the cover plate
5	Lift the cover plate away from the body of the instrument and remove it from the instrument cavity. Store it in a safe location away from the instrument. The reader is visible when the cover is removed.
6	<ul> <li>Open the door of the Plate Loading Tower and remove all of the plates from the tower.</li> <li>Note: Steps 6, 7, 8, 9, and 10 are the recommended best practice steps to remove and replace the lamp through the empty Plate Loading Tower. It is possible to remove and replace the lamp through the main body of the NEO Iris module. However the benefits of using your left hand through the Plate Loading Tower aperture are that: <ul> <li>a. You will reduce the risk of damage to the 4-probe pipettor arm.</li> <li>b. You will also be able to stand looking through the main body of the NEO Iris module so that you can visually coordinate the use of your left hand, and not just rely on the use of tactile hand navigation.</li> </ul> </li> </ul>

Step	Action
7	Using your left hand, reach into the top of the exposed camera reader module through the open and empty Plate Loading Tower and grip the exposed base of the camera lamp, just visible above the camera lamp socket.
	Note: The socket is to the right of the horizontal lamp installation assembly.         Pull the lamp horizontally to the left, away from the lamp socket.
	Note: The glass tip of the lamp on the left hand side of the lamp installation rests in an open ended bracket that suspends the bulb from the top of the camera module.

Step	Action
8	Once the malfunctioning lamp is detached from its socket, angle the bulb down to gently remove it from the open ended bracket. Remove the lamp from the instrument through the Plate Loading Tower.
	The brackets will now be empty (shown below).

Step	Action
9	Using your left hand, insert the new lamp through the Plate Loading tower and hook the glass tip of the lamp into the open ended bracket on the left hand side of the lamp installation.
	Make enough clearance on the right hand side so that the lamp base can be pushed horizontally into the lamp socket.
10	Push the lamp base firmly into the lamp socket.
	<b>Note:</b> the glass tip of the lamp should be held by the open ended bracket on the left hand side of the lamp installation.
11	Lift the cover plate onto the body of the instrument through the instrument cavity. The
	Internal cavity of the reader is now not visible when the protective cover is in place.
	slid into place at the base of the cover aperture.
12	Screw the two thumb screws clockwise to secure the protective plate in place.

Step	Action
13	Close the NEO Iris hood.
14	Power the NEO Iris on.
15	Contact Technical Support to determine whether additional activities are required.

## **Replacing the Y-pusher**

While performing error recovery after a plate jam, care should be taken when removing plates from the plate transport. If too much force is applied when removing the plates, it is possible that the Y- pusher may separate from the plate transport. Therefore it must be returned back to its correct location for correct NEO Iris operation to occur.

#### About the Y-pusher

The Y-pusher is located on the plate transport. The Y-pusher is the arm that moves the plates onto and off the transport.



#### Procedure

To replace the Y-pusher:

Step	Action
1	Power off the NEO Iris module.
2	Open the NEO Iris hood.
3	Move the transport to an accessible position.  Note: You must always use the silver knob on the top of the transport to
	move the transport around. Do not use any other part of the transport.
4	Position the Y-pusher with the teeth facing to the right and the 3 prongs facing downwards, so that the pusher looks like an upside down L.

Step	Action
5	Gently slide the Y-pusher fully back into position by pushing it with your fingertip.
	As you push, the Y-pusher should move freely, with some mechanical resistance, through the
6	Close the NEO Iris hood.
7	Power on the NEO Iris module.
	<b>Note:</b> After replacement of the Y-pusher, monitor the NEO Iris for transport related errors. Contact Technical Support if you need additional assistance.

# Chapter 11: Troubleshooting the NEO Iris

# In This Chapter

The NEO Iris is designed to automate error-free processing of blood samples. Nevertheless, errors can occur in practice. This chapter contains instructions on how to overcome such errors with minimum loss of assay results. Help screens are available in the software to aid you in error situations.

CHAPTER 11: TROUBLESHOOTING THE NEO IRIS	11-1
The Troubleshooting Process Steps	11-2
Using Error Codes to Troubleshoot	11-4
Troubleshooting Software Failure	11-10
Pipettor Self Check Failures	11-14
Clot Detection Recovery Process	11-18
Troubleshooting Plate Transport Errors	11-19
Troubleshooting Pipettor Errors	11-27
Troubleshooting Centrifuge Errors	11-32
Troubleshooting Incubator Errors	11-44
Troubleshooting Washer Errors	11-50
Troubleshooting Camera Reader Errors	11-58
Troubleshooting 14-lane and 5-lane Bay Errors	11-60
Troubleshooting Plate Tower Errors	11-62

# The Troubleshooting Process Steps

Troubleshooting consists of implementing a practical, systematic approach to problem solving. This approach focuses on:

- Observing, recognizing, and categorizing symptoms
- Identifying the probable cause(s)

The following troubleshooting model describes a four-step approach to defining symptoms, identifying problems, and implementing solutions. When troubleshooting the instrument you should also include considerations appropriate to your specific environment.

### **Identify the Problem**

The first step in the troubleshooting approach is to identify the error. It is important to clearly state the error or problem with cause unknown.

#### Data Gathering

The second step in troubleshooting is data gathering. The observations, data and clues about the problem are collected in this step, without attributing any one clue or observation as the cause.

Four questions must be asked:

- What? For example, what operation was being performed at the time of the error?
- *When*? For example, when did the error occur?
- *Where*? For example, where did the error occur?
- *Who***?** For example, who was operating the instrument when the error occurred?

It is important to make sure that data and observations be collected without attempting to determine the cause of the problem. All relevant data will be evaluated in the next step. It can be helpful to ask the following questions:

- ► What object is affected?
- ▶ What other objects could be affected but are not?
- ▶ Where do you see the problem?
- ▶ Where else could I expect to see this problem occur?
- ▶ Where on the object/instrument does the error occur?
- ▶ Where else on the object/instrument could this problem occur?
- When did the problem first occur?

- ▶ When else could the problem have occurred?
- ▶ When in the process flow does the problem occur?
- ▶ When else in the process flow has the error occurred?
- ► Is the problem repetitive?
- ► Is the problem random?
- ▶ In what pattern does the problem occur?
- ► How else could the pattern repeat itself?
- Does the problem occur for all operators?
- Does the problem occur for a specific operator?

#### **Evaluate Cause**

Once all relevant data, clues and observations have been gathered an evaluation can be performed. If you have experience with the problem it is acceptable to prioritize clues and observations to generate likely causes. Depending on the amount of data and the error being encountered it may be important to implement various evaluation methods.

Ask the following questions:

- Based on the data, what appears to be the likely cause?
- Based on my experience, what clues appear to be a likely cause?

#### **Action Plan**

Once the evaluation step is completed and a likely cause has been identified, the final step in troubleshooting an error is to implement an action plan. This step often consists of a standard operating procedure or some other step wise instruction that addresses or remedies the likely cause that was isolated in the previous step. Should the action plan fail to produce the expected change in the problem it would be necessary to return to the data gathering step and continue until the actual likely cause is identified and remedied

# Using Error Codes to Troubleshoot

Error codes display when the instrument needs to provide you with some additional information or if something is causing the instrument to function improperly. You can apply the previously described troubleshooting steps.

## Identify the Problem

The first step is to read the error message.

As a best practice, it is recommended that the operator press F12 on the keyboard to silence the alarm and then press Print Screen on the keyboard to print the error message for future reference. The printed error message can be used as a tool by Technical Support to help make an accurate diagnosis.

It is essential that you determine what the error is before trying to recover from it. For example, in the error message below, you can determine the following:



#### Data Gathering

Observe the error code and instrument to gather information about the failure. Ask guestions appropriate to the troubleshooting process to gather specific relevant data.

### **Evaluate Cause**

Based on observations and the data gathered during the previous step it is important to isolate the likely cause of the hardware failure. The error code generated above indicates a hardware failure. The last two to three words of the error code normally identify the step at which the error occurred. The recovery message may also give clues to the likely cause of the error.



Note: The recovery messages for some errors instruct the operator to "Call Service". If you are able to recover from the error then it is not necessary to call Technical Support.

#### **Action Plan**

When acting to fix a hardware failure you must first clear any jam or blockade before initializing the module in question to prevent further errors and complete error recovery. You must manually clear the error and complete the step at which the error occurred or let the instrument make another attempt at completing the step at which the error occurred before continuing work.



**Warning:** During error recovery, the process control allows you to open the safety hood to access components manually. You must be cautious during manual instrument deck interventions, because while you resolve an error in one system module, other operations may remain active.

After completing your assessment of the likely cause, determining an action plan, and clearing any blockade, Error Recovery must be performed. The following steps will guide you through the Error Recovery process.

Step	Action
1	Press the <b>OK</b> button in the error message dialog. This silences the audible alarm, if it has not already been silenced by using the <b>F12</b> key. The <i>Machine is Stopped!</i> dialog is displayed on the screen, allowing access to the <i>Error recovery</i> dialog, via the <b>Error recovery</b> button.

Step	Action
2	Always press the <b>Error recovery</b> button in the <i>Machine is stopped!</i> window to enter the <i>Error Recovery</i> window and, from there, perform manual or software supported error recovery actions.
	Description Transport - Y- Position Not Reached at position 08h params 01h while getting plate Recovery: Mask the plate cartier involved in the created Possible reasons: mechanical blockade in y-axis, y- adapting and wrong, y-encodersignals wrong Select module Reader Reader Contribute Washer1 Incubator1 Pipettor Loading Bay Tower
	The <i>Error Recovery</i> window provides initialization and recovery functions for each module. This window allows you to recover from the error or to abort the run immediately. <b>While</b> <b>you are in this window you will be able to open the safety hood without signaling the</b> <b>audible alarm and to manually manipulate the affected modules.</b> The <i>Description</i> and <i>Recovery</i> fields provide guidance as to the potential causes at the top
	of the window. The first word of the error message lists the module the error affected.

Step	Action
3	Press the appropriate module button from the <i>Select module</i> area. This will display recovery buttons and processes for the selected item.
	Decorption Pecovery: Mode the plate cases exolved in the case! Possible recorr: machanical blockade in y exis, y Select module FaccodeReader Reader Centraloge Internation Internatio Inte
	Remember that, in most cases, the first word in the error message tells you which module you need to work with. Refer to the print screen you obtained of the original error if necessary.
	You should clear any visible blockade or plate jam before initializing the module, or you may cause further errors. Further errors could delay the error recovery process and ultimately exceed processing limits for some assay steps, leading to invalidated plate results.
	Once you select the module you will be working with, several buttons will become visible on the right side of the error recovery screen.
	The number of recovery buttons and processes vary between modules and are used to initiate movement of the components, or initialize the module. Some modules may only allow you to initialize the module, while others may also allow you to move a module to a specific position or to open and close a door.
	Initializing the module should be the first button selected in this window, after clearing any blockades. Once initialization is complete, other options may become available, if necessary, to complete error recovery. During the error recovery process, the instrument should be returned to the step at which the error occurred or the step immediately following the step at which the error occurred. This will depend on what was manually performed and what actions were performed when the module was initialized.
	The individual module error recoveries are described later in this chapter.          Image: If the error occurred in the pipettor or washer module you should always select the Abort plate option after initializing the module in <i>Error recovery</i> .

Step	Action
4	When you are finished in the <i>Error recovery</i> dialog, press the <b>Close</b> button in the upper right corner of the window.
	Close

Step	Action
5	You will again see the <i>Machine is stopped!</i> dialog, but note that the two middle buttons are now available to you:
	<b>Try last step again</b> should be selected if you have not completed the step in which the error occurred either manually or via the software. This will signal the instrument to make another attempt to complete the step. For example, if the plate was being pulled onto the transport, and following initialization the plate is still not on the transport, <b>Try last step again</b> will signal the instrument to make another attempt to place the plate carrier onto the transport. Once the step is successfully completed, you will need to press <b>Continue work</b> . <b>Continue work</b> should be selected if:
	<ul> <li>You have manually completed the step in which the error occurred. For example, the transport experienced an error while getting the plate and, after resolving the error, you manually placed the plate onto the transport.</li> <li>You have already selected <b>Try last step again</b> and that step was successfully completed by the instrument.</li> <li>You have already selected <b>Try last step again</b> and no movement was made by the instrument.</li> </ul>
	<b>Abort run</b> should be selected if you wish to abandon all plates in process or if you cannot recover from the current error. If this is selected, you should reinitialize the instrument before scheduling more tests.
	If <b>Continue work</b> is selected, all assays will resume processing and the affected plate will be flagged with the error. For any plate, the results will be invalid if the mechanical error occurred at the pipettor or washer modules or the plate remains in the incubator for more than 60 minutes.

Step	Action
	Hemagglutination assay plate results will be invalid if the plate is prevented from going from the centrifuge to the reader within a software encoded time, because of the risk that negative results, having suspended red blood cells, will have cells settling out and therefore cause reader interpretation problems.
	Solid Phase Red Cell Adherence plate results will be invalid if, following the pipetting of indicator red cells, the plate is prevented from going into the centrifuge within a software encoded time, because of the risk of indicator red cells being neutralized and therefore cause the subsequent risk of false negative results.
	Capture plate results will be invalid if, following the plate wash step, the plate is prevented from going to the Pipettor station within a software encoded time, because of the risk of monolayer drying, which could cause erroneous results.

# Troubleshooting Software Failure

Although the instrument software is designed to minimize the occurrence of failures, you cannot entirely avoid software failures. The instrument design allows the hardware to continue working for some time, even when the main PC software is no longer running. This feature is dependent on the module activity at the time of the software failure.

An example of a module activity that will continue in the event of software failure is the incubation phase of the incubator. This feature provides a chance to recover from a software failure without losing the assay run during which the failure occurred.



**Note**: You can try to complete the current run after software failure, but you must reboot the PC and power cycle the instrument before new runs can be started.

# **Resolving Software Failure Errors**

To resolve software failure errors:

Step	Action	
1	If the software fails, the system displays one of two error messages, the <i>Program Error</i> window or the <i>Microsoft Visual C++ Debug Library</i> window.	
	Program Error  mpa exe has generated errors and will be closed by Windows. You will need to restart the program. An error log is being created.  Cancel	
	When the system displays the Program Error window, do not press Cancel. After the system	
	creates an error log, the <b>Cancel</b> button changes to <b>OK</b> . You may then press <b>OK</b> . The system	
	closes the message box and returns to the normal user interface.	
2	Wicrosoft Visual C++ Debug Library         View of the system of the system can cause an assertion failure, see the Visual C++ documentation on asserts.         (Press Retry to debug the application)         Abort       Retry         Ignore         Press the Retry button in the Microsoft Visual C++ Debug Library window. The system creates an error log. The system closes the message box and returns to the normal user interface.	
3	<ul> <li>The system displays a <i>Create Diagnostic Archive</i> icon on the Windows desktop. Double-click this icon to automatically generate an archive of all information available about this software failure. In addition, you must make a note of the following information:</li> <li>Date and time of software failure.</li> <li>Description of operator actions before the error occurred.</li> <li>Description of instrument actions before the error occurred.</li> <li>Call Technical Support with this information regarding the software failure.</li> </ul>	

Step	Action	
4	After collecting this information, restart the computer. This clears any partial software components still active in the memory and automatically restarts the main application. Because of the time these actions take, the instrument may have stopped moving and will automatically reactivate.	
5	After logging into the software, the message regarding <i>Data consistency</i> will appear indicating that the software was improperly closed. All runs in progress will complete (if the were not performing a critical step) and the system will require a full shutdown and initialization prior to new runs being started. Press the <b>OK</b> button.	
	Warning (Number: 00010105)         03.03.2010-08:34:23         Message:       Image: Imag	
6	However, because of the time delay, it is possible that the system will display a <b>Pipettor</b> <b>answer time-out</b> error message or there may be a visible stoppage in processing without an error message. To clear this message, you must press the <b>OK</b> button to enter the error message. If there is a visible stoppage, the instrument and software must be shutdown and initialized to continue.	
7	You must then immediately press the <b>Close</b> button to exit the Error recovery window.	

Step	Action		
8	Press the <b>Continue work</b> button on the <i>Machine is stopped!</i> dialog.		
	Failed step 03h - 10C8h		
	Error at module Transport		
	Error recovery		
	Try last step again		
	Continue work		
	Abort run		
	In many cases, the instrument continues the current run and provides results for all plates		
	can continue to occur, and you may be forced to abandon the run. If the instrument cannot		
	recover following a software failure, it either does not start when you press the <b>Continue</b>		
	work button, or stops again after a few actions.		
	Any runs that do complete after the software had been unexpectedly closed will have a warning flag associated with them indicating this. The results are valid; however it is recommended the plates are viewed and verified for accuracy.		
	Warning: After the current run is finished, you MUST shut down the software again, turn off the instrument, and then restart. If the attempt to complete the		
	run fails, you must also shut down both the software and the instrument, and		
	instrument and the PC.		

# Pipettor Self Check Failures

#### **Expected Results**

When the Pipettor Self Check maintenance task is performed and the pass status is achieved, there are wells that are expected to be flagged and also wells that are expected to not contain flags. The expected sequence is illustrated in the screen image below.



Column	Expected Flags
1	Not Enough Liquid or No Liquid flags are acceptable
2, 4, 5, 7	No flags should be present
3 and 6	Not Enough Liquid Flags only are acceptable – No Liquid flags will cause the task to fail

During the task, each instrument probe pipetting activity can be mapped to a collection of wells. The pipetting map is described in the table below.



**Note:** The number 4 left probe is the probe furthest away from the front of the instrument and number 1 is nearest to the front.

Mapping	Left Probes	Right Probe
Row A (columns 1-4)	4	
Row B (columns 1-4)	3	

Row C (columns 1-4)	2	
Row D (columns 1-4)	1	
Row D (columns 5-7)		Right probe

### **Troubleshooting Failures**

In the event of a failed status for the Pipettor Self Check task, you must view the plate flags to investigate the possible reasons for failure. There are many different reasons why the task can fail and from a troubleshooting perspective, the pattern or sequence of additional unexpectedly flagged or un-flagged wells can indicate the most likely causes of the failure. Some possible causes of failure are listed in the table below with some recommended recovery steps.

**Note:** The final step of any recommended recovery is to repeat the Pipettor Self Check task. You must contact Immucor Technical Support to report the problem if the recommended recovery is unsuccessful and the consequence is a persistently failing Pipettor Self Check task.

Column Number and Feature	Possible Reason for Failure	Recommended Recovery
No flag is generated for at least one of the wells in <b>column 1</b>	Aspiration failure – Possibly due to a blocked probe or failed syringe.	Inspect the syringes for leakage and replace as needed. Perform Initialization and any additional required maintenance as prescribed.
Flags are present in any of the wells of <b>column 2</b>	Air in the system liquid lines.	Initialize the instrument. Verify the connections of the pump rack location (below the syringes) if the air cannot be removed by priming. Reconnect any loose connections and repeat the system liquid priming.
	Syringe is damaged or is leaking and it may not be delivering enough liquid to the wells. Probes taught too high and unable to aspirate all of the	Inspect the syringes for leakage and replace as needed. Perform Initialization and any additional required maintenance as prescribed. Verify Check Pipettor Reference positions are acceptable and verify the PCheck tool was

Column Number		
and Feature	Possible Reason for Failure	Recommended Recovery
	liquid in the tool.	properly loaded in transport frame
No flag is generated for at least one of the wells in <b>column 3</b>	Aspiration failure – Possibly due to a blocked probe or failed syringe.	Inspect the syringes for leakage and replace as needed. Perform Initialization and any additional required maintenance as prescribed.
No Liquid Flag generated in <b>column</b> <b>3</b>	Probe(s) taught too high and unable to detect the remaining liquid	Verify Check Pipettor Reference positions are acceptable and verify the PCheck tool was properly loaded in transport frame
Flags are present in any of the wells of <b>column 4</b>	Air in the system liquid lines (which might be the result of a leak from the pump rack connections or the de-bubbler mechanism).	Initialize the instrument. Inspect the connections of the pump rack location (below the syringes) and also the de-bubbler mechanism if the air cannot be removed by priming. Reconnect any loose connections and repeat the system liquid priming.
	A rinse pump has failed and is not dispensing system liquid.	Contact Immucor Technical Support to report the problem.
	Too much foam is generated during the rinsing step.	Contact Immucor Technical Support to report the problem.
Flag is present in well <b>D</b> of <b>column 5</b>	Air in the system liquid lines.	Initialize the instrument. Verify the connections of the pump rack location (below the syringes) if the air cannot be removed by priming. Reconnect any loose connections and repeat the system liquid priming.
	Syringe is damaged or is leaking and it may not be delivering enough liquid to the wells.	Inspect the syringes for leakage and replace as needed. Perform Initialization and any additional required maintenance as prescribed.
No flag is generated for well <b>D</b> in <b>column</b> 6	Aspiration failure – Possibly due to a blocked probe or failed syringe.	Inspect the syringes for leakage and replace as needed. Perform Initialization and any additional required maintenance as
Column Number and Feature	Possible Reason for Failure	Recommended Recovery
--	--	--
		prescribed.
No Liquid Flag generated in <b>column</b> <b>6</b>	Probe(s) taught too high and unable to detect the remaining liquid	Verify Check Pipettor Reference positions are acceptable and verify the PCheck tool was properly loaded in transport frame
Flag is present in well <b>D</b> of <b>column 7</b>	Air in the system liquid lines (which might be the result of a leak from the pump rack connections or the de-bubbler mechanism).	Initialize the instrument. Inspect the connections of the pump rack location (below the syringes) and also the de-bubbler mechanism if the air cannot be removed by priming. Reconnect any loose connections and repeat the system liquid priming.
	A rinse pump has failed and is not dispensing system liquid.	Contact Immucor Technical Support to report the problem.
	Too much foam is generated during the rinsing step.	Contact Immucor Technical Support to report the problem.

# Clot Detection Recovery Process

If a probe detects a clot during sample aspiration, the system displays an error message alerting you to the event, and the probe remains stationary above the tube in question. This is your opportunity to inspect the tip of the probe for any hanging clots. If visible hanging clots need to be removed, lift up the NEO Iris hood while the error message is still being displayed and use an absorbent dry wipe to remove the hanging clot and clean the tip of the probe. If no hanging clots are visible, use either the **Skip** or **Continue** button to proceed.

Error (Number: 00020D28) 23.03.2010-13:16:17	
Message: COP: Clot detected on the left pipettor arm, processing is halted. Affected Probe [Probe #, Liquid ID -> Target Well]: 1, F694 b -> D4	
Recovery Message:           Remove clot manually.           "Skip" will not dispense with the affected probes.           "Continue" will process the affected probes.	)
Skip Continue	

You can press either the **Skip** or **Continue** button to clear the error message and cause the NEO Iris to proceed according to two (2) different protocols.

Skip	Continue
------	----------

By pressing the **Skip** button, the probe is sent to the wash tower for cleaning before processing the remaining on-board samples. The sample that had contained the clot is bypassed and no results will be generated for that sample.

By pressing the **Continue** button, the probe will retry the sampling of the tube that was identifed as containing the clot. Repeating this aspiration step may cause another clot detection error. The clotted sample may be taken off the NEO Iris for further testing by an alternative methodology if required.

# Troubleshooting Plate Transport Errors

### Introduction

This section describes errors in the plate transport. They include

- Symptom 1: The plate transport does not get a plate out of a module.
- Symptom 2: The plate transport fails to put a plate into a module.



**Note:** If a plate transport error occurs, then it is recommended to mark the plate carrier involved with an identifying mark so that if additional plate transport errors occur in the future, the initial plate carrier can be tracked to determine if the causative factor is related to a specific defective plate carrier.

## **Common Causes of Plate Transport Errors**

You may encounter various issues that can prevent the plate transport from working properly. The causes are separated into two types: those you may encounter during routine operation of the instrument, and those you may encounter during the error recovery process. You must be aware of these causes and the impact they can have on your NEO Iris.

Common causes during normal operation include, but are not limited to:

- Failure to snap the plate completely into the carrier (crashes at the reader)
- Strips that are not snapped completely into the frame (possible crash at reader or incubator)
- Object or material dropped into the path of the transport (impedes transport movement)
- Use of broken plate carriers (loose/missing pieces can impede movement of transport or Y-pusher)
- Failure to place the plate completely into the plate loading tower without closing the door (Y-pusher cannot reach plate to remove from tower)
- Loading a plate into the plate loading tower while the transport is accessing the tower (plate moves too far and impedes the transport as it moves in the X direction)

Common causes for additional errors during error recovery include, but are not limited to:

- Applying too much pressure to the Y-pusher when removing plates (causing misalignment of the Y-pusher)
- Not clearing a plate jam before initializing the transport (transport crashes when it begins to move)
- Not selecting **Init transport** when initializing the transport (transport may not be in the correct position to **Continue** or **Try last step again**)

### Transport Fails to Remove a Plate from a Module

### **Error Message**

Transport: Y-position not reached at position 08h params while getting plate

I	Error (Number: 14030106) 15.09.2009-10:42:01	
Message:		
Transport: Y- Position Not Re while getting plate	eached at position OBh params O1h	
Recovery Message: Mark the plate carrier involve mechanical blockade in y-axi encodersignals wrong	d in the crash! Possible reasons: s, y-motorsignals wrong, y-	
		×

### Potential Root Cause(s)

- Too much friction from the plate position in the module
- Transport frame broken or twisted
- Transport teach position at this module is sub-optimal
- Y-pusher may need to be cleaned

#### Procedure

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.

Step	Action
3	The software is displaying the <i>Error Recovery</i> window. Press the <b>Transport</b> button.
	This will display recovery buttons and processes for the selected item.
	Transport recovery         Init transport         1. Identify what motion failed. Manually complete or undo the failed motion.         2. Press 'Init transport'. Wait until the small motions in the leftmost position have finished.         3. Close error recovery. Press 'Try last step again'.         3.1 If a 'Plate Position Error' occurs return to 'Error recovery.' and 'Transport', press 'Init transport' again. Wait long enough! Close error recovery.         4. Check that the transport stopped moving. Then press 'Continue work'.         5. After the run, mark the plate carrier that was involved in the error. If this plate carrier recurrently causes errors replace it.
	Note:The Init Y-Pusher button can be used to initialize the Y-pusher. TheDodge right and Dodge left buttons can be used to slightly adjust the position
	of the transport. The Service position button can be used to move the position
	of the transport to the serviceable position in front of the incubator module.

Step	Action		
4	Determine where the transport module is currently located on the instrument (Step 1 of on- screen instructions)		
	<b>Scenario 1:</b> The plate carrier is still completely inside the position that it should have been taken from and the Y-pusher is partially or fully engaged in the plate carrier <u>OR</u> the plate carrier is partially pulled out of the position. The Y-pusher has encountered resistance during its motion.		
	a. If the transport is at the reader, incubator, or either pipettor stations, open the hood to better access the transport		
	b. <u>Carefully</u> grab the <i>Silver Transport Knob</i> (to the right of where the plate sits on the transport) and pull upwards slightly		
	c. Push the plate fully onto the module or station that the transport is trying to pick it up		
	d. Manually retract the Y-pusher back into the transport module		
	e. Verify the transport can move freely		
	f. Press the Init transport button (Step 2 of the onscreen instructions) and then follow the remaining on-screen <i>Transport recovery</i> instructions. Instructions are also available in step 3 of this procedure.           Image: Silver Knob         Si		
	<b>Note:</b> If the plate transport must be moved, you must only hold the silver knob on top of the transport to move its position.		

Step	Action
	<b>Scenario 2:</b> The plate carrier is still completely inside the position that it should have been taken from. The Y-pusher is not engaged in the plate carrier, because it encountered resistance while moving toward the plate carrier. The plate carrier is fully pulled out of its position and located almost on the plate transport. The carrier encountered resistance in the last few millimeters of the action.
	a. Press the <b>Init transport</b> button (Step 2 of the onscreen instructions) and then follow the remaining on-screen <i>Transport recovery</i> instructions.
	<b>Scenario 3:</b> The plate carrier is fully pulled out of its position and located almost on the plate transport. The carrier encountered resistance in the last few millimeters of the action.
	a. Press the <b>Init transport</b> button (Step 2 of the onscreen instructions).
	Do not press the Try last step again button because the failed step (i.e. unloading
	a. the plate onto the transport) has already been performed.
	Close <i>Error recovery</i> and then press the <b>Continue work</b> button.



**Note:** Refer to **Chapter 10 – Maintaining the NEO Iris** for information regarding the Y-pusher.

# Transport Fails to Load a Plate into a Module

#### **Error Message**

Transport: Y-position not reached at position 08h params while getting plate

Em	03.03.2010-08:50:33	oj
Message:		
Transport: Y-Position Not Reach while releasing plate	ed at position 07h params 0	1h 🔺
łecovery Message: Mark the plate carrier involved in mechanical blockade in y-axis, y-	the crash! Possible reasons: motorsignals wrong, y-	
encodersignals wrong		y z
	ОК	

### Potential Root Cause(s)

- Plate not seating fully in transport frame
- Transport frame broken or twisted
- Too much friction from the plate position in the module
- Transport teach position at this module is sub-optimal
- Y-pusher may need to be cleaned

### Procedure

Step	Action	
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then Press the <b>OK</b> button on the message window.	
2	The system displays the <i>Machine is stopped!</i> dialog. Press the Error recovery button.	
3	The software is displaying the <i>Error recovery</i> window. Press the <b>Transport</b> button.         This will display recovery buttons and processes for the selected item.         Init vectors       1 identify what motion failed. Manually complete or undo the failed motion.         Init vectors       2 Press Init transport. Warually complete or undo the failed motion.         2 Press Init transport.       2 Press Init transport.         3 Doce error recovery return to Error recovery e	

Step		Action
4	Detern screen	nine where the transport module is currently located on the instrument (Step 1 of on- instructions)
	<b>Scena</b> resista	rio 1: The plate carrier is still almost completely on the transport system. It encountered nce while moving toward the module position
	a.	Press the <b>Init transport</b> button (Step 2 of the onscreen instructions) and then follow the remaining on-screen <i>Transport recovery</i> instructions. Instructions are also available in step 3 of this procedure.
	<b>Scena</b> is part the mo	<b>rio 2:</b> The plate carrier is halfway inside the destination module position. The Y-pusher ially or fully engaged in the plate carrier <u>OR</u> The plate carrier is almost completely in odule position.
	a.	If the transport is at the reader, incubator, or either pipettor stations, open the hood to better access the transport.
	b.	<u>Carefully</u> grab the <i>Silver Transport Knob</i> (to the right of where the plate sits on the transport) and pull upwards slightly.
	c.	Push the plate fully onto the module or station that the transport is trying to pick it up.
	d.	Manually retract the Y-pusher back into the transport module.
	e.	Verify the transport can move freely.
	f.	Press the Init transport button (Step 2 of the onscreen instructions).
	g.	Do not press the <b>Try last step again</b> button because the failed step (i.e. loading the plate onto the destination module) has already been performed.
	h.	Close <i>Error recovery</i> and then press the <b>Continue work</b> button.
	F	<b>Note:</b> If the plate transport must be moved, you must only hold the silver knob

on top of the transport to move its position.

Step	Action
	Scenario 3: The Y-pusher is no longer engaged with the plate carrier. It encountered resistance in the last few millimeters of the action
	a. Verify the plate is pushed fully onto the destination module
	Press the <b>Init transport</b> button (Step 2 of the onscreen instructions) and then follow the remaining on-screen <i>Transport recovery</i> instructions. Instructions are also available in step 3 of this procedure.



**Note:** Refer to **Chapter 10 – Maintaining the NEO Iris** for information regarding the Y-pusher.

# Troubleshooting Pipettor Errors

# Introduction

Various types of pipettor errors can occur, involving either the XYZ arm movements or the pumps. However, the recovery process to continue the current run is always the same. If the error is persistent, you are required to perform different actions to fix the various types of problems.

The complexity of the pipetting sequences renders it impossible for the system's error recovery routines to know exactly at which substep of a move the error occurred. Therefore, no **Try last step again** function is available for the pipettor, and no valid results can be expected from a plate that suffered a pipettor error.

The NEO Iris has two pipetting stations, left and right. Therefore two plates are frequently being pipetted simultaneously. For example, samples are transferred to a plate on the left side at the same time as reagents are dispensed onto a second plate on the right side. After a pipettor error, both arms need to be re-initialized, so the processing of both of these plates is disrupted.

After re-initializing the pipettor, it is recommended that you abort the processing of both affected plates. You can do this with the **Abort plate** button on the error recovery screen (if two plates are present, both are aborted). If you select this button, the remainder of the current pipetting sequence is finished, but then the system returns the plate(s) to the loading tower.



Warning: You should always select the Abort plate button (not the Abort run button).This will abort the plate with the error as well as the plate disrupted by re-initialization.Failure to select the Abort plate button allows processing to continue for a plate that may not be valid. If you select the Continue work button without aborting the plate, the system

will invalidate all results generated from the plate that suffered the pipettor error.

### **Pipettor Error Message**

The *Message* field informs you

- That the pipettor is involved.
- Of the drive or component that failed.
- Of the error mode (open loop) during which the error occurred.

The *Recovery Message* field lists some possible reasons for failure and instructs you to call Technical Support for service.

The system can generate other error messages when the pipettor fails. The potential messages are:

Error Message	Description
Left (or Right) X (or Y or Z) – Position not Reached.	The pipettor did not reach the position to which it was supposed to move. A mechanical blockage probably impeded the movement of the pipettor in the respective direction.
Left (or Right) Arm Move would cause Collision.	If the respective arm moves, the two arms would collide. Call Technical Support if this error occurs.
Move to undefined target.	The pipettor is sent to a target that is unknown to the pipettor. Call Technical Support if this error occurs.
Motor Error – Open Loop (Left or Right) X (or Y or Z) Overload.	The pipettor did not reach the position to which it was supposed to move. A mechanical blockage probably impeded the movement of the pipettor in the respective direction.
Diluter Error.	An error from the diluter pump. Call Technical Support if this error occurs.
Left (or Right) diluter answer timeout!	The dispense pumps did not communicate with the main processor on the instrument for an extended period of time.
	To rectify this error, open the <b>Error recovery</b> dialog and press the <b>Close</b> button (without initializing the pipettor), then press the <b>Continue work</b> button. The NEO Iris should proceed without causing an error.

## Procedure

The initial recovery attempt to continue the current run is always the same, except for two error messages, *Diluter answer timeout* and *Left (or Right) Diluter answer timeout*. These can be recovered without causing a flag. All of the other causes require different solutions to resolve the underlying problem.

The table below describes how to troubleshoot a pipettor error:

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then Press the <b>OK</b> button on the message window.

Step	Action	
2	The system displays the <i>Machine is stopped!</i> dialog.	
	Press the <b>Error recovery</b> button.	
3	The hood can be opened if needed. The system displays the <i>Error Recovery</i> window. Press the <b>Pipettor</b> button. This will display recovery buttons and processes for the selected item.           Pipettor recovery         Check       1. Press 'Check reference' to check for proper probe alignment.         2. Press 'Abort plate'.       3. Error recovery closes automatically and processing commences.         Remark: Never use 'Try last step again' after a pipettor error.	
4	Check the NEO Iris to see where the pipettor probes are located. Lift the probes carefully upward out of any tubes or off of metallic surfaces. If the probes have visible red cells on them, gently clean each probe with an alcohol pad prior to continuing. Check to see if a probe was bent. If so, bend it back into a straight position at the very top of the probe while supporting the black probe adaptor with one hand.	
5	Follow the on-screen <i>Pipettor recovery</i> instructions.	
6	Press the <b>Check reference</b> button to verify the probe positions via the check pipettor reference maintenance task.	
	Note: Refer to Chapter 10 – Maintaining the NEO Iris for more information regarding how to check the pipettor reference.	
7	The <b>Abort plate</b> button becomes activated after the pipettor has been reinitialized.	
	Press the <b>Abort plate</b> button to abandon processing of the plates that are on the pipetting stations.	

Step	Action	
8	A pop-up dialog indicating the plate(s) that aborted will appear and an audible alarm will be generated. The plate(s) will automatically be returned to the plate loading tower during the next scheduled transport step.	
	<b>Note:</b> There may be a delay before the plate is returned to the tower because the NEO Iris will wait until the time allotted for pipetting has expired.	

# Troubleshooting Centrifuge Errors

## Introduction

Centrifuge recovery requires several manual steps, because centrifuge loading, centrifugation, and unloading consist of several substeps, each of which can fail. The **Try last step again** and **Continue work** buttons function on the substep level.

Centrifugation consists of the following system substeps and will be referenced below in the specific error recovery steps:

- 1. Place the plate on the load unit (a transport system step).
- 2. Open the centrifuge door.
- 3. Push the plate carrier into the centrifuge.
- 4. Close the centrifuge door.
- 5. Perform one or several centrifugation processes or shaking steps and brake afterwards.
- 6. Open the centrifuge door.
- 7. Pull the plate carrier from the centrifuge.
- 8. Close the centrifuge door.
- 9. Take the plate from the load unit (a transport system step).

For a successful recovery, the centrifuge, door, and load table subcomponents must ALL be in the appropriate states.

## Centrifuge Fails to Initialize During Complete Initialization

#### **Error Message**

#### Centrifuge: Rotor Drive Error

### Potential Root Cause(s)

- The main power for the centrifuge is OFF
- Something is obstructing the rotor movement
- Mechanical failure
- Failed recovery from previous mechanical error

#### Procedure

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.
2	Press the <b>Abort run</b> button.
3	Cancel the initialization process by pressing the <b>Cancel</b> button in the current initialization window.
4	Verify that the mains power switch for the centrifuge is turned <b>ON</b> (the switch must be depressed in the I position). Note: The centrifuge power switch is located at the rear of the centrifuge near the power cord. Centrifuge power switch

Step	Action	
5	Begin the NEO Iris initialization process by pressing the initialization main menu bar button and confirming the initialization process on the resulting confirmation dialog.	
6	Troubleshooting is complete if the initialization completes without error.	
	If the error occurs again, continue to step number 7.	
7	Turn the centrifuge power <b>OFF</b> (the switch must be depressed in the <b>O</b> position).	
	Note: The centrifuge power switch is located at the rear of the centrifuge, next to the centrifuge power cord.	
8	Remove the Centrifuge Loading Unit (CLU) cover (see item <b>A</b> in photo below), if it has not already been removed. Use a 2.5 mm Allen key to remove the Allen head screw (and its associated washer) that secures the CLU cover.	
	A: CLU cover.	
	<b>B</b> : Centrifuge service hatch cover.	
	a. Use a wrench to remove the four (4) huts securing the service hatch cover (see item <b>B</b> in photo above).	
	b. Remove the service hatch cover by lifting it away from the four permanent screw posts. You can now access the inside of the centrifuge through the hatch.	

Step	Action
9	If there is a plate in the centrifuge, remove the white plate frame only from the rotor. The transport carrier cannot be removed at this time. If there is no plate in the centrifuge, skip to step number 10.
10	Push the transport frame as far into the centrifuge (away from you) as possible to ensure it is not the cause of the error.
11	Manually rotate the centrifuge rotor completely around to ensure there are no obstructions. Call Technical Support if the rotor will not rotate completely and the obstruction cannot be identified.
12	Install the service hatch.
	a. Hang the service hatch cover back onto the four permanent screws.
	<ul> <li>Replace the four nuts securing the hatch cover by hand tightening them over the four associated screw posts and then use a wrench to securely tighten the nuts (see item <b>B</b> in photo above).</li> </ul>
13	Install the CLU cover and use a 2.5 mm Allen key to replace the Allen head screw (and its associated washer) that secures the CLU cover.
14	Turn the centrifuge power back <b>ON</b> (the switch must be depressed in the I position).
15	Initialize the NEO Iris by pressing the initialization main menu bar button and confirming the initialization process on the resulting confirmation dialog.
16	Contact Technical Support if the error persists.

# Plate Carrier Is Not Loaded Properly into Centrifuge

### **Error Messages**

Centrifuge: Door Movement Error

Centrifuge: Plate present after loading

#### Potential Root Cause(s)

- Alignment issue
- Something is obstructing the plate loading
- Mechanical failure

If the plate fails to load, the plate position will dictate the error generated. If the plate fails early in the process, the Plate Present after loading message will be generated. If the plate fails towards the end of the process, the Door movement error will be generated because the plate is beyond the sensor and the door is trying to close. Usually a considerable amount of mechanical noise accompanies this error.

### Procedure

Step	Action	
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.	
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.	

Step	Action	
3	The system displays the Error Recovery dialog. Press the Centrifuge button.	
	This will display recovery buttons and processes for the selected item.	
	Centrifuge recovery         Init centrifuge         1. Take note of the following system properties:         1.1 Where is the microplate, inside the centrifuge or on the loading table?         1.2 Is the centrifuge's door opened or closed?         Transfer         position         3. Press 'Init centrifuge'.         0pen door         4. Return the system to its previous state.         4.1 Put the plate on the table. Only if it was fully inside the centrifuge put it there.         Close door         5. Close error recovery, press 'Try last step again'.         5.1 After the motion has finished, press 'Continue work'.         5.2 If 'Try last step again' leads to another error, then repeat error recovery and press 'Continue work'.         Unload plate	
	Text is displayed in the <i>Centrifuge recovery</i> area describing symptoms to look for in the diagnosis of the problem. Follow the on-screen <i>Centrifuge recovery</i> instructions.	
4	Remove the Centrifuge Loading Unit (CLU) cover. Use a 2.5 mm Allen key to remove the Allen head screw (and its associated washer) that secures the CLU cover. Remove the plate carrier, even if it is already partially or completely inside the centrifuge.	
5	Press the <b>Init centrifuge</b> button. The CLU and centrifuge initialize, including attempting to unload a plate from the centrifuge. Wait until this is completely finished.	
6	Return the plate carrier to the loading table.	
7	Press the <b>Open door</b> button and then press the <b>Load plate</b> button,	
	Contact Technical Support if this step fails.	
8	Press <b>Close</b> to leave error recovery. The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Continue work</b> button. Press the <b>OK</b> button on the verification dialog to indicate that affected modules have been initialized. The door closes. Centrifugation and test processing continues and the system closes the <i>Machine is stopped!</i> dialog after a few seconds.	

## Plate Carrier Is Not Unloaded Properly from Centrifuge

#### **Error Message**

### Centrifuge: No Plate present after loading

### Potential Root Cause(s)

- Alignment issue causing gripper to lose plate in process
- Something is obstructing the plate unloading
- Mechanical failure

#### Procedure

Step	Action		
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.		
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.		
3	The system displays the <i>Error Recovery</i> dialog. Press the <b>Centrifuge</b> button. This will display recovery buttons and processes for the selected item. Centrifuge recovery Internifuge 1. Take note of the following system properties: 1.1 Where is the microplate, inside the centrifuge or on the loading table? 1.2 Is the centrifuge's door opened or closed? 2. Remove the microplate, press Transfer position' if necessary. 3. Press Init centrifuge'. 4. Return the system to its previous state. 4. Peturn the system to its previous state. 5. Close error recovery, press Try last step again'. 5.1 Mart the molon has finished, press "Continue work'. 5.2 If Try last step again' leads to another error, then repeat error recovery and press "Continue work'. Text is displayed in the <i>Centrifuge recovery</i> area describing symptoms to look for in the diagnosis of the problem.		
	Follow the on-screen Centrifuge recovery instructions.		

Step	Action	
4	Remove the Centrifuge Loading Unit (CLU) cover. Use a 2.5 mm Allen key to remove the Allen head screw (and its associated washer) that secures the CLU cover.	
	Remove the plate carrier, even if it is partially or completely inside the centrifuge.	
	Note: If the plate has difficulty sliding out of the centrifuge, press the <b>Transfer</b> position button to attempt a re-alignment of the load position. Repeat this attempt to manually remove the plate.	
5	Press the <b>Init centrifuge</b> button. The load unit and centrifuge initialize, including attempting to unload a plate from the centrifuge. Wait until this is completely finished.	
6	Return the plate carrier to the loading table.	
7	Press the <b>Open door</b> button.	
8	Press the <b>Close</b> button to leave error recovery. The system displays the <i>Machine is stopped!</i>	
	dialog. Press Continue work. Press the OK button on the verification dialog to indicate that	
	affected modules have been initialized. The door closes. Centrifugation and test processing	
	continues and the system closes the <i>Machine is stopped!</i> dialog after a few seconds.	

# **Centrifuge Stalls**

The system can generate other error messages when centrifugation stalls. The potential messages are:

Error Message	Description
Centrifuge: Error while moving vibration adjust motor.	Before actual centrifugation, movement of imbalance adjustment is impossible, there is friction, or there is a motor problem.
Centrifuge: Vibration adjust movement out of range or internal error.	Before or during centrifugation, movement of imbalance adjustment reached limits, plate weight is out of range, or centrifuge adjustment has changed.
Centrifuge: Vibration adjust movement lasts too long.	During centrifugation, the balance point cannot be found, the centrifuge is not mounted correctly, or mechanical play causes imbalance.

Error Message	Description
Centrifuge: Imbalance sensor value after imbalance compensation still too high.	During centrifugation, the balance point is found, but vibration level is too high, the centrifuge is not mounted correctly, or mechanical play causes imbalance.
Centrifuge: Vibration adjust sensor detected too much imbalance.	During centrifugation, the vibration level is above allowed limit, the centrifuge is not mounted correctly, or mechanical play causes imbalance.
Centrifuge: Rotor drive error: Open loop.	At the beginning of centrifugation: A possible motor or sensor failure. During centrifugation: A probable mechanical blockage or the plate carrier moved out of table.
Final centrifugation velocity not reached.	There is partial motor failure, power driver failure, or friction in bearings.

Although all of these different causes may require different solutions to resolve the underlying problem, the initial recovery to try and continue the current run is always the same.

### Procedure

Step	Action	
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.	
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.	

Step	Action
3	The system displays the Error Recovery dialog. Press the Centrifuge button.
	This will display recovery buttons and processes for the selected item.
	Centrifuge recovery         Init centrifuge         1. Take note of the following system properties:         1.1 Where is the microplate, inside the centrifuge or on the loading table?         1.2 Is the centrifuge's door opened or closed?         2. Remove the microplate, press 'Transfer position' if necessary.         3. Press 'Init centrifuge'.         4. Return the system to its previous state.         4.1 Put the plate on the table. Only if it was fully inside the centrifuge put it there.         Close door         5. Close error recovery, press 'Try last step again'.         5.1 After the motion has finished, press 'Continue work'.         5.2 If 'Try last step again' leads to another error, then repeat error recovery and press 'Continue work'.         Unload plate         Text is displayed in the <i>Centrifuge recovery</i> area describing symptoms to look for in the
	diagnosis of the problem.
	Follow the on-screen <i>Centrifuge recovery</i> instructions.
4	Verify that the mains power switch for the centrifuge is turned <b>ON</b> (the switch must be depressed in the I position). <b>Note:</b> The centrifuge power switch is located at the rear of the centrifuge, next to the centrifuge power cord.
	Centrifuge
5	Press the Init centrifuge button.
	If the centrifuge initializes without error, skip to step number <b>17</b> .
	If the error persists, continue to step number <b>6</b> .
6	Turn the centrifuge power <b>OFF</b> (the switch must be depressed in the <b>O</b> position).

Step	Action
7	The Centrifuge Loading Unit (CLU) cover (see item <b>A</b> in photo below) should already be removed. If it is not, remove it now.
	Use a 2.5 mm Allen key to remove the Allen head screw (and its associated washer) that secures the CLU cover.
	A: CLU cover.         B: Centrifuge service hatch cover.
	a. Use a wrench to remove the four (4) nuts securing the service hatch cover (see item <b>B</b> in photo above).
	b. Remove the service hatch cover by lifting it away from the four permanent screw posts. You can now access the inside of the centrifuge through the hatch.
8	Only remove the white plate frame from the rotor. The transport carrier cannot be removed at this time.
9	Push the transport frame as far into the centrifuge (away from you) as possible to ensure it is not the cause of the error.
10	Manually rotate the centrifuge rotor completely around to ensure there are no obstructions. Call Technical Support if the rotor will not rotate completely and the obstruction cannot be identified.

Step	Action
11	Replace the service hatch.
	a. Hang the service hatch cover back onto the four permanent screws.
	<ul> <li>b. Replace the four nuts securing the hatch cover by hand tightening them over the four associated screw posts and then use a wrench to securely tighten the nuts (see item <b>B</b> in photo above).</li> </ul>
12	Turn the centrifuge power back ON (the switch must be depressed in the I position).
	• -
13	Press the Init centrifuge button.
	If the centrifuge initializes without error, continue to step number 14.
	Contact Technical Support if the error continues.
14	Place the plate on the CLU if it was previously removed.
15	Install the CLU cover and use a 2.5 mm Allen key to replace the Allen head screw (and its associated washer) that secures the CLU cover.
16	Press the <b>Open door</b> button and then press the <b>Load plate</b> button. Press the <b>Close door</b> button.
17	Press the <b>Close</b> button to close the <i>Error recovery</i> dialog. The software displays the <i>Machine is stopped!</i> dialog.
	Press the <b>Try last step again</b> button. The centrifuge will attempt to perform the centrifuge step that generated the initial error.
18	Wait until the centrifugation or shake sequence is complete and then press the <b>Continue work</b> button.
	Press the <b>OK</b> button on the dialog to confirm that the affected module(s) have been initialized. The remaining assay sequences will be performed.

# Troubleshooting Incubator Errors

### **Incubator Door Does Not Close**

#### **Error Message**

Incubator: Flap jam while closing (open/close light barrier error)

#### Potential Root Cause(s)

- There is something obstructing the incubator door from closing
- The sensor detecting a closed door state is defective

#### Procedure

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then Press the <b>OK</b> button on the message window.
2	The system displays the <i>Machine is stopped!</i> dialog.
	Press the <b>Error recovery</b> button.

Step	Action
3	The system displays the <i>Error recovery</i> window. Press the <b>Incubator1</b> button.
	Select module       Reader       Reader       Centrifuge         Washer1       Incubator1       Pipettor         Transport       Sample       Tower         This will display recovery buttons and processes for the selected item.
	Incubator recovery     Service       Init incubator     1. If necessary remove obstacle.       2     Press 'Init incubator'.       Stot No.     3. If applicable use slot open/close to achieve required state.       Open slot     4. Close error recovery. Press 'Try last step again'. then 'Continue work'.       Close slot     Close slot
4	Check to see whether anything is obviously jamming the door. The plate carrier might not be fully inserted into the incubator slot or friction in the slide mechanism of the door might be blocking its motion. Remove any visible obstruction. Push the plate fully into the incubator.           Note:         It may be necessary to disengage the transport Y-pusher from the plate frame if that is blocking the door movement.
5	Press the Init incubator button. All incubator bay doors should close during initialization.
	Contact Technical Support if this step fails.
6	Press the <b>Continue work</b> button.
	Press the <b>OK</b> button on the dialog to confirm that the affected module(s) have been initialized. Test processing continues and the system closes the <i>Machine is stopped!</i> dialog after a few seconds.

### Incubator Door Does Not Open

#### **Error Message**

Incubator: Flap jam while open (open/close light barrier error)

#### Potential Root Cause(s)

- There is something obstructing the incubator door from opening
- The door is sticking and the friction was too great
- The sensor detecting a closed door state is defective

#### Procedure

Step	Action	
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then Press the <b>OK</b> button on the message window.	
2	The system displays the <i>Machine is stopped!</i> dialog. Press Error recovery.	
3	The system displays the <i>Error Recovery</i> window. Press the <b>Incubator1</b> button. This will display recovery buttons and processes for the selected item.	
4	Press the <b>Init incubator</b> button. All incubator bay doors should close during initialization. Contact Technical Support if this step fails.	
5	Press the <b>Close</b> button to leave error recovery. The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Try last step again</b> button. The incubator door that previously failed to open should open at this time.	

Step	Action	
6	Press the <b>Continue work</b> button.	
	Press the <b>OK</b> button on the dialog to confirm that the affected module(s) have been initialized. Test processing continues and the system closes the <i>Machine is stopped</i> ! window after a few seconds.	

## Wrong Incubator Door Opens

This results in a transport system error because the unexpectedly closed incubator door is a physical blockage to the transport.

• This can occur when two incubator doors are stuck together or if the door above the desired position falls down

#### Procedure

Step	Action	
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then Press the <b>OK</b> button on the message window.	
2	The system displays the <i>Machine is stopped!</i> dialog. Press Error recovery.	
3	The system displays the <i>Error recovery</i> window. Press the <b>Transport</b> button. This will display recovery buttons and processes for the selected item.           Init ransport         1. Identify what motion failed.           Init transport         1. Identify what motion failed.           Init transport         1. Identify what motion failed.           Init transport         2. Press Thin transport.           Init Y-Pusher         3. Obse error recovery.           Press Try last step again.         3.1 If a Plate Position Error occurs           Service position         3.1 If a Plate Position Error occurs           Dodge left         0. Check that the transport stopped moving.           Dodge left         0. After the run, mark the plate carrier that was involved in the error.           It this plate carrier recurrently causes errors replace it.	
4	Determine and take note of the slot that the transport was trying to access. It is most likely to be the one just under the open slot. It is possible that the flap directly above the opened slot	
5	Press the <b>Init transport</b> button to initialize the transport. The transport slowly moves to the bottom left initialization position.	

Step	Action
6	While the transport is initializing, press the Incubator1 button. This will display recovery buttons and processes for the selected item.         Incubator recovery       Service         Initiacubator       Press 'Init incubator'.         Stot No       If applicable use slot open/close to achieve required state.         Open slot       Close error recovery.         Press 'Try last step again'.         then 'Continue work'.
7	<b>Optional:</b> Manually slide the incorrectly opened door up to open the correct slot. If it does not stay in place, further action is necessary. Sticky liquid at the connection edge between the doors can cause this dropping of the door above the desired one. Clean the edges between the doors with absorbent material.
8	Enter the slot number ( <b>Slot No</b> .) for the incubator bay that needs to be opened. Note: Doors are numbered 1 (bottom-most) – 15 (top-most).
9	Press the <b>Close</b> button. The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Try last step again</b> button.
10	Press the <b>Continue work</b> button. Press the <b>OK</b> button on the dialog to confirm that the affected module(s) have been initialized. The software closes the <i>Machine is stopped!</i> dialog after a few seconds.

# Troubleshooting Washer Errors

The complexity of the washer sequences renders it impossible for the system's error recovery routines to know exactly at which strip position of a move the error occurred. Therefore, the **Try last step again** function should not be used for the washer, and no valid results can be expected from a plate that suffered a washer error.

# Liquid Handling Errors

The potential washer liquid handling messages and the corresponding troubleshooting steps are:

### Liquid Overflow

#### Error Message

#### Washer1: Liquid overflow Sensor in overflow trough detected liquid

Error (Number: 2C03020 26.02.2010-14:19:44	)B)
Message:	
Washer1: Liquid overflow Sensor in overflow trough detected liquid	*
	⊻ <b>(</b> )
Recovery Message:	
Check for aspiration blockage. Initialise washer.	
	y.
DK	

#### Potential Root Cause(s)

- Blocked aspiration channel
- Poor manifold positioning
- Strips not in expected positions for maintenance assays

#### Procedure

Step	Action	
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.	
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.	
3	The system displays the <i>Error recovery</i> window. If the transport is in front of the washer, press the <b>Transport</b> button and then press the <b>Service position</b> button.	
	<b>Note:</b> This will move the transport to a position in front of the incubator and out of the way of the washer.	
	Do not perform this step if the transport is not in front of the washer.	
4	Press the <b>Washer1</b> button. This will display recovery buttons and processes for the selected item.	
	Washer recovery         Init washer         2. Press 'Init washer'.         Abort plate         3. Press 'Abort plate'.         4. Error recovery closes automatically and processing commences.	
5	Press the <b>Init washer</b> button (step 2 of the on screen recovery instructions). The manifold slowly moves to the left side of the washer.	
6	Manually remove the plate from the washer module.	
7	<b>Completely</b> dry the overflow trough using an absorbent towel. If an additional liquid overflow error occurs, initialize the washer again and continue drying the trough.	
8	Return the plate to the washer module in the correct position	
9	Press the <b>Abort plate</b> button.	
	The software will automatically continue processing. The aborted plate will be returned to the tower by the transport module during the next scheduled transport step.	
10	If it has been determined that the overflow error is a result of a blocked manifold channel, the manifold should be removed and cleaned according to the instructions in <b>Chapter 10</b> - <b>Maintaining the NEO Iris</b> .	

### **Dispense Check Failure**

#### Error Message

Washer1: Dispense check failure:

Top Probes ???????

#### Bottom Probes ???????

Error (Number: 2C03) 26.02.2010-14:18:13	D20F)
Message:	
Washer1: Dispense check failure : 1=failed row HGFEDCBA Top Probes 10000001 Bottom Probes 10000001	
Recovery Message: Please check for proper priming, wash head and dispens Initialise washer, abort plate. Report to support.	
ОК	

#### Potential Root Cause(s)

- Blocked dispense channel
- Poor seal between manifold and washer assembly
- Air bubbles in the tubing

#### Procedure

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.
Step	Action
------	---
2	The system displays the <i>Machine is stopped!</i> dialog. Press the Error recovery button.
3	Use the Print Screen to analyze the error.
	a. The values displayed for the Top Probes and Bottom Probes should be analyzed.
	b. A value of '0' indicates the expected liquid has been detected in the well(s) and no error condition exists.
	c. A value of '1' indicates the expected liquid was not detected in the well(s) and an error condition does exist. These are the wells that must be investigated.
	d. In the example above, the expected volume was not detected in channels A and H.
4	Press the Washer1 button. This will display recovery buttons and processes for the selected item.          Item.         Init washer         Press 'Init washer'.         Abort plate         Press 'Init washer'.         Press 'Init washer'.         Abort plate         Error recovery closes automatically and processing commences.
5	Press the <b>Init washer</b> button (step 2 of the on screen recovery instructions). The manifold slowly moves to the left side of the washer.
6	Press the <b>Abort plate</b> button.
	The software will automatically continue processing. The aborted plate will be returned to the tower by the transport module during the next scheduled transport step.
7	If it has been determined that the dispense error is a result of a blocked manifold channel, the manifold should be removed and cleaned according to the instructions in <b>Chapter 10</b> - <b>Maintaining the NEO Iris</b> .

#### **Aspiration Check Failure**

#### Error Message

Washer1: Aspiration check failure:

1=failed row HGFEDCBA

Top Probes ???????

#### Bottom Probes ???????

Error (Number: 2C0302 26.02.2010-14:13:39	20E)
Message:	
Washer1: Aspiration check failure: 1=failed row HGFEDCBA Top Probes 10000000 Bottom Probes 10000000	
Recovery Message: Please check wash head and vacuum pump Initialise washer, abort plate. Report to support.	
ок	

#### Potential Root Cause(s)

- Blocked aspiration channel
- Poor seal between manifold and washer assembly
- Vacuum pump failure
- Check valve failure

#### Procedure

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.
3	Use the Print Screen to analyze the error.
	a. The values displayed for the Top Probes and Bottom Probes should be analyzed.
	b. A value of '0' indicates that all expected has been aspirated from the well(s) and no error condition exists.
	c. A value of '1' indicates that all expected has <u>not</u> been aspirated from the well(s) and an error condition does exist. These are the wells that must be investigated.
	d. In the example above, the too much liquid remained after aspiration and was detected in channel H.
4	Press the Washer1 button. This will display recovery buttons and processes for the selected item.
5	Press the <b>Init washer</b> button (step 2 of the on screen recovery instructions). The manifold slowly moves to the left side of the washer.
6	Press the <b>Abort plate</b> button. The software will automatically continue processing. The aborted plate will be returned to the tower by the transport module during the next scheduled transport step.
7	If it has been determined that the aspiration error is a result of a blocked manifold channel, the manifold should be removed and cleaned according to the instructions in <b>Chapter 10</b> - <b>Maintaining the NEO Iris</b> .

### **Mechanical Problems**

The system can generate other error messages when mechanical problems occur. The potential messages are:

Error Message	Description
Washer 1: Command () microplate missing.	There is no plate in the transport frame or the frame sensor is defective. Verify a microplate is loaded and call Technical Support if initialization does not resolve issue.
PC: Washer – horizontal (x-motor) positioning error.	The motor that moves the manifold over the plate could not reach the intended position. Verify whether anything, such as tubing, is hindering manifold movement.
Washer1: X movement error lost 17 steps. Mechanical blockage or drive damaged	If no physical blockage is obvious and the error is sporadic, contact Technical Support.
PC: Washer – vertical (z-motor) positioning error possibly incorrect height of bottom touch is detected:	The upward or downward movement is impeded by something like a mechanical blockage. Verify whether the manifold is in the correct position.
	If this error occurs at the priming trough, visually verify that the manifold aspiration/dispense needles are not hitting the side of the trough. This may be due to a displacement of the trough.
	Verify all of the test strips are seated in the plate frame and the plate frame is fully pressed into the transport frame.

Although the listed root causes can require different solutions to resolve the underlying problem, initial recovery is always the same to try to continue the current run.

### Procedure

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.
3	The system displays the <i>Error recovery</i> window. Press the <b>Washer1</b> button. This will display recovery buttons and processes for the selected item.          Washer recovery         Init washer         2. Press 'Init washer'.         3. Press 'Abort plate'.         4. Error recovery closes automatically and processing commences.
4	Press the <b>Init washer</b> button (step 2 of the on screen recovery instructions). The manifold slowly moves to the left side of the washer.
5	Press the <b>Abort plate</b> button. The software will automatically continue processing. The aborted plate will be returned to the tower by the transport module during the next scheduled transport step.

# Troubleshooting Camera Reader Errors

### Introduction



**Note:** Camera reader recovery is a special case, because the transport processor controls the reader door. Therefore, the error module indicated is Transport. However the *Reader recovery* button must be utilized.

Errors in the image reading process itself, such as wrongly positioned plates or no light, generate an error message and no useful images. The affected plate will automatically be aborted. The machine does not enter the *Stopped* state. No recovery procedure is required. You need to investigate persistent image reading failures and eliminate their root cause.

If you must stop assay processing because the image reading persistently fails, it is recommended that the current schedule is **Canceled** to abort the processing. This leads to abandoning the assay processing completely while preserving unused reagent volume.

### **Reader Door Does Not Open or Close**

#### **Error Message**

#### Transport: Reader Door Error (open or close)

#### Potential Root Cause(s)

- Movement of door obstructed
- Reader door motor has failed
- Reader door sensors have failed

#### Procedure

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.

Step	Action	
3	Press the <b>Reader</b> button. This will display recovery buttons and processes for the selected item.	
4	Press the <b>Init reader</b> button.	
5	Press the <b>Open door</b> and <b>Close door</b> buttons a few times to see whether the door is moving freely. If a jam is evident, try to remove the cause.	
6	Press the <b>Open door</b> or <b>Close door</b> buttons, depending on the initial error.	
7	Press the <b>Close</b> button to close the error recovery dialog. The system displays the <i>Machine is stopped</i> ! dialog. Press the <b>Try last step again</b> button.	
	Press the <b>Continue work</b> button. Press the <b>Yes</b> button on the dialog to confirm that the affected module(s) have been initialized. The reader loads or unloads the plate carrier and the door closes. The system closes the <i>Machine is stopped!</i> dialog after a few seconds.	

# Troubleshooting 14-lane and 5-lane Bay Errors

### Introduction

There are some errors that can occur in the 14-lane and 5-lane bays. The software displays these errors in the *Event Log* (most loading bay errors do not result in pop-up dialogs). These errors are mainly related to the barcode scanning process.

### Barcodes on Rack Could Not Be Interpreted

The rack appears with question marks and the Event Log error message informs you:

- That either the 14-lane or 5-lane bay and the PC are involved.
- That a barcode on the rack could not be interpreted.

### **Position Barcode Missing**

If a rack is loaded and there is an issue reading one of the positional barcodes on the rack, the rack loading will be denied and a message will appear in a pop-up dialog and an entry will be made to the *Event Log*.



Note: The message will not appear in the log list scroll seen from the main screen.

- Remove the rack and inspect the positional barcodes
- If error repeats with the same rack, attempt to load a different rack

Warning (Number: 0E010 03.03.2010-10:25:49	043)	
Message:		
Sample Loading Bay: Position barcode missing	Ā	
Recovery Message:		
Inspect rack barcode labels	Ă	Ŷ
ОК		

#### Procedure

Step	Action	
1	Follow the steps for <i>Rescanning Unread Sample Barcodes</i> published in <b>Chapter 6</b> – <b>Instrument Testing Operation</b> .	
	<b>Note:</b> When reinserting the rack, you must make sure that your fingers do not visually obstruct barcodes, including the rack ID barcode.	
2	If the reader still cannot read the barcodes, remove the rack again and check to see if any of the barcodes on the rack are removed, dirty, incomplete, or otherwise unreadable. Replace the rack if necessary.	
3	Close the loading bay dialog and allow the scanner to initialize. A message will appear in lower log list.	
4	Open the loading bay dialog and attempt to scan the rack again. If the issue persists, initialize the instrument and retry the scan.	
5	Contact Technical Support if the issue cannot be resolved.	

# Troubleshooting Plate Tower Errors

The plate tower errors relate to the opening and closing of the plate loading back shutter. The shutter is a mechanism which helps to consistently manually load plates into the tower so that the transport can use the same coordinates to remove plates. Prior to the transport loading or unloading from the tower, the shutter will be opened. Operators should not load plates during this time.

### Plate Loading Shutter Will Not Close

#### **Error Message**

#### Tower: Plate Loading Shutter will not close

#### Potential Root Cause(s)

- Plate pushed into tower too far by operator while transport was accessing the module
- Transport did not push the plate far enough into the tower when loading
- Movement obstructed or mechanical failure

#### Procedure

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then Press the <b>OK</b> button on the message window.
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.
3	Open the tower door if it is currently closed.

Step	Action
Step 4	Action         The system displays the Error Recovery dialog. Press the Tower button.         Select module         BarcodeReader Reader Centrifuge         BarcodeReader Reader Centrifuge Tower       Pipettor         Washer1 Incubator1 Pipettor       Difference         Transport Loading Bay Tower       Tower         This will display recovery buttons and processes for the selected item. Follow the on screen instructions.         Tower recovery       1. Open tower front door.         Open shutter       1. Open tower front door.         2. Pull all plates out by 1cm.
	Llose shutter       3. Press 'Close shutter'.         4. Push all plates softly against shutter.         5. Close error recovery. Press 'Try last step again', then 'Continue work'.

### Plate Loading Shutter Will Not Open

#### **Error Message**

#### Tower: Plate Loading Shutter will not open

#### Potential Root Cause(s)

- Plate pushed into tower too far by operator while loading plates
- Movement obstructed or mechanical failure

#### Procedure

Step	Action
1	Press the <b>F12</b> key to stop the alarm. Obtain a <b>Print Screen</b> of the message and then press the <b>OK</b> button on the message window.
2	The system displays the <i>Machine is stopped!</i> dialog. Press the <b>Error recovery</b> button.

Step	Action
3	Open the tower door if it is currently closed.
4	The system displays the <i>Error Recovery</i> dialog. Press the <b>Tower</b> button.          Select module       Reader         BarcodeReader       Reader         Incubator1       Pipettor         Transport       Sample         Loading Bay       Tower         This will display recovery buttons and processes for the selected item. Follow the on screen instructions.         Tower recovery       1. Open tower front door.         Qpen shutter       2. Pull all plates out by 1cm.         Close shutter       3. Press 'Close shutter'.         4. Puch all plates softly against shutter.       5. Close error recovery.         Press Try last step again', then 'Continue work'.       Press Try last step again', then 'Continue work'.

# Chapter 12: Limitations of Use and Warnings

# In This Chapter

This chapter is a collected list of all of the limitations of use and warnings that are located within the body of the text where they are most relevant to the information written. Each chapter and appendix is sequentially listed below with the limitations of use and warnings assigned to it.

CHAPTER 12:	LIMITATIONS OF USE AND WARNINGS	12-1
Limitations	of Use	12-2
Warnings		12-16

# Limitations of Use

### Limitations by Chapter

The table below lists the Limitations of Use of the NEO Iris that are contained in the following chapters:

- Chapter 1 Introduction to the NEO Iris
- Chapter 2 Hardware Components
- Chapter 3 System Software Navigation
- Chapter 4 Security
- Chapter 5 Instrument Start-Up
- Chapter 6 Instrument Testing Operation
- Chapter 10 Maintaining the NEO Iris
- Appendix A Preparing the NEO Iris for First Use
- Attachment 1 for NEO Iris Operator Manual

# Chapter 1 – Introduction to the NEO Iris

lcon	Description
IVD	Limitation: The Neo is for in vitro diagnostic use.
	Limitation: Barcodes can be no longer than 18 characters in length.

lcon	Description
	<b>Limitation</b> : Barcodes must have a module size larger than 0.2mm and a bar width ratio between 2.25:1 to 3:1.
	Decodability grade (grading system of A to F; A being best, F is failing) of C or better is required for consistent reading of barcodes on the instrument. This grade measures the bar width consistency throughout the barcode label. It is usually an indication of print quality of the barcode label.
	The minimum number of characters in the barcode is 3.
	The length of the barcode is variable but must be completely visible with a quiet zone (white space on each end of the label) of 3 mm when placed in the sample or donor rack.
	The minimum height of the barcode must be 10 mm.
	If barcodes have parameters outside of these specifications, barcode misreads can occur on the instrument.
	Pipe characters () are not permitted as part of a sample identification.
	<b>Limitation</b> : If samples have barcode identification information that is eighteen (18) characters in length and the first three (3) characters are identical to the first three (3) characters of the assay control material in assays which include plate or run controls, then the sample will be interpreted as a replicate of the control material. In this case, either the plate will fail unnecessarily (if the sample reacts differently than expected for the corresponding control), or the plate will pass but no results for that sample will be produced (if the sample reacts as expected for the corresponding control). Such a condition can also be exhibited when an assay such as crossmatch or antigen screening assay generates a set of circumstances such that the combination of donor and primary sample barcode identification begins with the same first three (3) characters as one of the control material barcode identifications.

	<b>Limitation</b> : A site visit by an Immucor representative is required to configure the Z position (downward) on the instrument for the C racks and the specific small-volume pediatric tubes in use at your site. C racks cannot be used on the instrument without this configuration. Sample probe crashes will occur without this configuration, when used in conjunction with the C racks. If differently sized small-volume pediatric tubes are subsequently used after the Z position configuration is performed using the originally designated small-volume pediatric tubes, re-configuration may be required to prevent possible sample probe crashes
	into the bottom of the new tubes.
	<b>Limitation</b> : Time stamps for instrument activity may not be accurate around Daylight Saving Time (DST) when a given activity spans a time period falling on both sides of the actual change of time for DST. The following recommendations are published to provide guidance on how to mitigate these time stamp inaccuracies. Allow assays to finish if they are already running during the DST change and remove the racks once processing is complete. Do not interact with the instrument (e.g. loading plates or starting assays) during the DST time change. Initialize the instrument after the DST time change is finished prior to beginning any further assays.

# **Chapter 2 – Hardware Components**

lcon	Description
	<b>Limitation</b> : Immucor requires the use of phosphate buffered (approximately 15mM) isotonic saline, pH 6.5-7.5 (PBS), on the NEO Iris system. Reactions between an antibody and its antigen may be weakened if acidic or unbuffered saline is used.
	Using saline and/or deionized water in PBS preparation from sources with systems in place to control proliferation of microbes helps to reduce the chance for microbial bioburden on the system. Excessive microbial bioburden can cause degradation of system or assay performance.
	<b>Limitation</b> : At least 250 $\mu$ l of packed red blood cells need to be present in a sample tube to ensure that the probe picks up red blood cells and not plasma (only for those assays that require red blood cells).
	At least 500 $\mu$ l of plasma or serum needs to be present in a sample tube to ensure that the probe picks up plasma or serum.
	<b>Limitation</b> : Laboratory ambient temperature and humidity affects the room temperature incubation bays, and an elevated ambient room temperature may disrupt assays that need to be incubated at specified temperature ranges, as published in the relevant package insert.
	<b>Limitation</b> : The NEO Iris must be switched on at least 30 minutes prior to the first plate read to allow the reader lamp to warm up. Reading of plates prior to completion of this warming period can cause incorrect negative reading of weakly positive reactions.

# **Chapter 3 – System Software Navigation**

lcon	Description
	<b>Limitation</b> : The timestamps for assay plate events that occur simultaneously with the end of Daylight Saving Time (DST) may not accurately reflect the actual time of those events, and it is not possible to predict what these timestamps will display due to their erratic nature. However, in this instance, the assay is successfully completed and the results are unaffected. These erratic timestamps are limited to the end of DST. By running the instrument initialization after the end of DST, the time stamping process is corrected. The timestamps for assay plate events that occur simultaneously with the beginning of DST are accurate and do not exhibit the same erratic features as the timestamps at the end of DST.
	Limitation: Validate Method is otherwise known as reflex testing. Reflex testing can only be ordered once per assay in response to results for a given sample number. The instrument WILL NOT order reflex testing again if that sample (with the same barcode) is repeated with the same assay on the instrument at a later phlebotomy date. This may cause incomplete test results to be indicated as final and made available for export. For example, when a sample is run and assigned the result of A Pending, the reflex testing procedure orders the reflex assay of Weak_D. If the sample tests positive for Weak D, then the final result of A Weak D will be released. A few days later, another sample with the same barcode is tested on the instrument and the result of A Pending is assigned again. However, on this occasion, the instrument WILL NOT order the reflex assay of Weak_D, and this second result of A Pending will be the final result. The historic and latest final results are not identical.
	This limitation also applies to the reflexive ordering of antibody ID testing for antibody screen positive samples. The consequence is that the second antibody identification test in the history of a given sample number will not be ordered for the second antibody screen positive sample, and so on for subsequent samples. After the instrument database is archived and the data is removed from the database, the reflex testing procedure will again order a reflex assay on a given sample ID that requires it because the current sample ID is again new to the instrument software
	<u>Limitation</u> : After the reflex testing procedure begins on a sample ID, the value of the result element not being resolved is not updated. For example, if a sample ID is run using the

lcon	Description
	ReflexABO assay and is assigned the result of B NTD, the reflex testing procedure will order a repeat of the ReflexABO assay to resolve result element 2 (the Rh result). If the result from this repeated ReflexABO assay is O Rh Positive, the reflex testing procedure will continue and the current result would be B Rh Positive. The ABO result of O that was determined during the repeat testing is not considered for incorporation into the final result.
	<b>Limitation</b> : In the event of an instrument computer hard drive crash, all data will be lost from the computer. Any data that was not previously archived will be lost.
	<b>Limitation</b> : Other than the Test per Assay Report, the outputs from the <i>Statistics</i> tab within the Utilities dialog must be verified on site by the end user. Immucor makes no claims as to the accuracy of those outputs.
	<b>Limitation</b> : The camera reader will verify the presence of clean and correctly positioned strips on a plate before starting an assay. If the reader detects the presence of strips that should have been deselected, such as wrong position of used strips, then the plate is aborted.
	<b>Limitation</b> : If an incorrect, but valid, plate expiration date is entered and an assay using this plate lot is run, then the incorrect expiration date is stored for that lot in the software. If the incorrect date is then corrected in the software following this assay run, and this correct date is further into the future than the original incorrect date, the data management system cannot update the incorrect stored date. Consequently, if an assay is run using this plate lot again, at a date of testing between the original incorrect date and the corrected date, then these results will be displayed in the Plate View window with an associated warning message of <b>Plate failed criteria</b> , <b>please check Interpretation Failures</b> . The validity of these results is not compromised by this chain of events that generated this warning message, and these results can still be exported or printed.
	<b><u>Limitation</u></b> : Stopping the system with the <b>Emergency Stop</b> button interrupts instrument sampling.

# Chapter 4 – Security

The following Limitations of Use are contained in this chapter:

lcon	Description
	Limitation: It is not possible to:
	• Assign to a new user higher user rights than the current operator.
	• Edit the user rights or change the password of a user with higher user rights than the
	current user.
	• Add a user who is already present on the user list.

# Chapter 5 – Instrument Start-Up

lcon	Description
	<b>Limitation</b> : The instrument must be switched on at least thirty (30) minutes prior to the first plate read to allow the reader lamp to warm up. Reading of plates prior to completion of this warming period can cause incorrect negative reading of weakly positive reactions.
	<b><u>Limitation</u></b> : The <b>Use DMS Archive</b> check box control option allows you to view archived results that are saved on a disk using the DMS. You can only view archived results through the DMS by inserting the disk into the drive of the computer.
	<b>Limitation</b> : The COP: Serial Buffer Deleted error message may appear during initialization. In this case click the <b>OK</b> button. However, if this error occurs in other situations it is recommended to complete the current run(s) but do not start any new plates. Initialize the instrument using the Initialize button before starting new plates.
	<b><u>Limitation</u></b> : Once the Auto Logoff has occurred, only users with equal or higher access rights than the user that was logged in when the Auto Logoff was triggered are able to access the system.

# **Chapter 6 – Instrument Testing Operation**

lcon	Description
	<b>Limitation</b> : When a sample is actively being tested by an assay and the same assay is requested again for that same sample from the LIS, a second request is added. This will cause the sample to appear in the Resource Overview window for a duplicate test. To avoid this, open and close the Loading Bay dialog (via Load Samples of the Start Run Assistant) AFTER the LIS worklist orders have been received but BEFORE the Resource Overview window is opened.
	<b><u>Limitation</u></b> : Excluding QC assays, some assays will not be processed under the following conditions:
	A corQC EXTEND Standard vial loaded in the 5-lane bay will be recognized; however this reagent will be listed as missing in the <i>Resource Overview</i> window due to the risk of a probe arm crash.
	When running a DAT assay, a DAT Positive Control Cell vial loaded in the 5-lane bay will be recognized; however this reagent will be listed as missing in the Resource Overview window due to the risk of a probe arm crash.
	Under both of these scenarios, the assays will not be able to be scheduled until the necessary reagent is loaded into the 14-lane bay.
	<b>Limitation</b> : When you insert a rack into the activated lane, the system scans the barcodes. The system displays the barcode data and interpretation (reagent name, lot number, expiration date) on the right-hand side of the on-screen dialog. The system also displays the volume remaining in each container. In the case of Immucor lot numbered reagent vials, the volume is based on the full volume of the container (for new vials) or the remaining volume (if the vial was previously used on the instrument). If the actual volume is less than that displayed (for example, if the reagent vial has also been used for manual testing), you may enter a reduced volume. It is not possible to enter an increased volume for the bottle.

lcon	Description
	<b>Limitation</b> : If a reagent expires during a testing process, the results will not be flagged. For example, a reagent that will expire after midnight is loaded onto the instrument and scheduled to an assay before midnight of the expiration date. However, the test will not be completed until after midnight when that vial has expired. The instrument prevents any test runs with this vial in the future.
	<b>Limitation</b> : When running the Crossmatch assay, if there are not sufficient wells remaining on a plate to accommodate the scheduled number of crossmatch tests for a given sample such as a sample to be crossmatched with three donor units and only one well remains on a plate, then the three tests will be batched on the next plate, leaving the remaining well on the first plate unused. If subsequent samples are scheduled for crossmatch assay, then the instrument cannot backtrack to the unused well(s) on the previous plate such as the one unused well in this example. Under these circumstances, such wells will remain unused.
	<ul> <li><u>Limitation</u>: If the common waste container is full and the waste icon turns red during sample processing, you will lose the interrupted testing run only if critical processing points exceed acceptable time limits because of excessive time removing the waste.</li> <li>Examples of exceeding critical acceptable limits are:</li> <li>Too much time has elapsed between ABORh plate shaking and reading.</li> <li>A Capture-R Ready-Screen plate is incubated beyond the programmed time limit.</li> </ul>

# Chapter 10 – Maintaining the NEO Iris

lcon	Description
	<b><u>Limitation</u></b> : Maintenance actions for the NEO Iris verify that specific modules of the instrument are functioning at the required specifications. The actions described are critical to the NEO Iris assay performance.
	If a maintenance item is not successfully performed for a module within the required time interval, then any assays that use that module will not run. For example, all assays use the reader. If you do not perform the daily <b>Pipettor Self Check</b> task then all assays are locked down and no assays can run.
	<b><u>Limitation</u></b> : Before using any cleaning or decontamination methods except those recommended by the manufacturer, users should check with the manufacturer that the proposed method will not damage the equipment.
	<b><u>Limitation</u></b> : Immucor requires the use of phosphate buffered (approximately 15mM) isotonic saline, pH 6.5-7.5 (PBS), on the NEO Iris system. Reactions between an antibody and its antigen may be weakened if acidic or unbuffered saline is used.
	Using saline and/or deionized water in PBS preparation from sources with systems in place to control proliferation of microbes helps to reduce the chance for microbial bioburden on the system. Excessive microbial bioburden can cause degradation of system or assay performance.
	<b>Limitation</b> : The instrument must be switched on for at least thirty (30) minutes prior to take flatfield images so that the reader lamp can warm up.
	<b>Limitation</b> : The instrument must be switched on for at least thirty (30) minutes prior to beginning <b>Take Flatfield Images</b> maintenance task so that the reader lamp can warm up.

# Appendix A – Preparing the NEO Iris for First Use

<b><u>imitation</u></b> : Adequate space must be provided for ventilation and also accessing the instrument, both during operation and in the event of servicing. The cabinet has two (2) spacing pins which are 4.8 inches (12.1 cm) in length. The instrument must be positioned such that an overall clearance distance of 30 inches (76.2 cm) must be maintained between he back of the instrument and the wall. This 30 inches (76.2 cm) includes the length of the spacing pins.
mproper connection of mains supply
<ul> <li>mproper connection of the instrument and the peripheral devices to the mains</li> <li>supply can cause serious personal injury with potentially deadly consequences and material damage (e.g. fire).</li> <li>Only use grounded connection and extension cables with sufficient capacity (voltage and current) to connect the instrument and any peripheral devices to the mains power supply.</li> <li>Never remove ground connections.</li> <li>Grounding of the instrument and its peripheral devices to the same protective earth potential shall be ensured.</li> <li>The use of a multi-outlet power strip is not allowed!</li> <li>Only use power cables that fulfill the minimum requirements for this instrument.</li> </ul>
<b>imitation</b> : <b>Service</b> : The instrument should be serviced by Immucor authorized service bersonnel. Only qualified technical personnel should perform troubleshooting and service brocedures on internal components.
<ul> <li><u>cimitation</u>: Environmental conditions: Do not expose the instrument to temperature extremes. For proper operation, ambient temperature range should remain between 18°C and 30°C. Performance may be adversely affected if temperatures fluctuate above or below his range.</li> <li>Shipping and storage temperature range: 0°C to 40°C.</li> <li>Relative humidity range: 20% to 80% (non-condensing).</li> </ul>

lcon	Description
	<b>Limitation:</b> Before using any chemicals (e.g. decontamintion solutions) except those recommended by the manufacturer, users should check with the manufacturer that chemicals will not damage the equipment or create potential hazards.
	Limitation: Warranty: Failure to follow preventive maintenance protocols may void the warranty. Refer to Chapter 10 – Maintaining the NEO Iris.
	<b>Limitation</b> : Changes or modifications to this unit not expressly approved by the manufacturer could void the user's authority to operate the equipment.
	<b>Limitation</b> : Do not install additional software on the PC. This will void your warranty and service contract. Third-party software can interfere with the controlling software and result in a loss of sample data.
	<b>Limitation</b> : <b>Warranty</b> : You must save all of the original packaging materials. If you need to ship the instrument to Immucor for repair or replacement, the original packing must be used. Other forms of commercially available packing are not recommended and can void your warranty.
	<b>Limitation</b> : You must decontaminate the instrument before it is repackaged. Follow the procedures in <b>Chapter 10 – Maintaining the NEO Iris</b> to decontaminate the instrument. The responsible laboratory authorities must state on the shipping documents that the instrument has been decontaminated and provide a signature for verification. Immucor reserves the right to refuse to accept an instrument that does not have such a decontamination statement.

# **Attachment 1 for NEO Iris Operator Manual**

lcon	Description
	<b>Limitation</b> : At least 250 µl of packed red blood cells need to be present in a sample tube to ensure that the probe picks up red blood cells and not plasma (only for those assays that require red blood cells).
	the probe picks up plasma or serum.
	<b>Limitation</b> : Red blood cell samples only collected with EDTA anticoagulants can be tested on the instrument. The use of other anticoagulants for testing on the instrument must be verified in the corresponding reagent Instructions for Use. Serum samples can also be tested on the instrument for tests that do not require red cells. Samples obtained from tubes containing neutral gel separators may produce falsely positive results and should therefore not be tested on the instrument.
	<b>Limitation</b> : The dead volume, or the depth to which probes cannot reach, is 200 $\mu$ l of sample.
	<b>Limitation</b> : When loading a test tube, containing red blood cell donor unit segment contents, into a Donor rack for compatibility testing with the IgG_XM assay, the contents of the test tube must not exceed a liquid height of 1.3 inches, otherwise the clot detection mechanism can be falsely alarmed.

lcon	Description
	Limitation: Samples that exhibit excessive hemolysis should not be tested on the instrument. Samples that exhibit a hemolysis grade of 3+ or greater must not be tested on the instrument, because they may generate erroneous results. Refer to the photograph below for a hemolysis grade of 3+ (color guide). For assays using Capture-R® Select, do not use hemolyzed samples for creating a monolayer. Fragmented red blood cell membranes will interfere with monolayer formation.
	<b>Limitation:</b> Samples that exhibit excessive lipemia should not be tested on the instrument. Samples that exhibit a dosage above 600 mg/dL of lipids must not be tested on the instrument because they may generate erroneous results.
	<b>Limitation</b> :Samples that exhibit excessive icterus should not be tested on the instrument. Samples that exhibit a dosage above 30 mg/dL of bilirubin must not be tested on the instrument because they may generate erroneous results.
	<b>Limitation</b> : Instrument studies have demonstrated that the probe washes are sufficient to prevent carry-over of samples that have an antibody titer of up to 1024 (typical of those samples encountered in a blood bank setting). It is important to note that these studies apply only to routine patient or donor samples. Sample material that is provided in surveys has been shown to occasionally cause carry-over at lower titers (a result of the manufacturing process for these materials).
	<b>Limitation</b> : For assays using Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> and Capture-R <sup>®</sup> Ready-ID <sup>®</sup> microplates, avoid testing previously frozen specimens that have undergone multiple freeze-thaw cycles, as these may give erroneous or inconsistent test results.

#### Incidents related to the device:

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

# Warnings

### Warnings by Chapter

The table below lists the warnings for the Neo that are contained in the following chapters:

- Chapter 1 Introduction to the NEO Iris
- Chapter 2 Hardware Components
- Chapter 3 System Software Navigation
- Chapter 6 Instrument Testing Operation
- Chapter 7 Test Results
- Chapter 10 Maintaining the NEO Iris
- Chapter 11 Troubleshooting the NEO Iris
- Appendix A Preparing the NEO Iris for First Use

### Chapter 1 – Introduction to the NEO Iris

lcon	Description
	<b>Warning</b> : Blood samples, liquid waste, used microplates, and consumed liquid reagent containers contain potentially biohazardous material.
	<u>Warning</u> : Always wear protective gloves and clothing when handling blood samples, liquid waste, used microplates, or consumed liquid reagent containers. All blood samples, liquid waste, used microplates, and consumed liquid reagent containers must be discarded following the standard practice of the laboratory.
	<b>Warning</b> : All blood products must be treated as potentially infectious. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.
4	<b>Warning</b> : Follow basic electrical hazard awareness to reduce the risk of injury due to prohibited electricity exposure.

<u>Warning</u> : Follow all of the necessary precautions to prevent exposure to and potential injury from instrument mechanical movement. Keep all instrument protective covers in place when operating the instrument to reduce the risk of operator injury due to instrument mechanical movement.
<u>Warning</u> : Follow all of the necessary precautions to prevent exposure to and potential injury from barcode laser scanners. Do not look directly into the laser beam of scanners or any reflections of the beam from a mirror-like surface. Exposure to the laser beam light can cause eye damage and permanent injury.

# Chapter 2 – Hardware Components

lcon	Description
	Warning: Inadvertent operator collision with the cabinet doors or the pull out cabinet shelf can cause operator injury.
	<b>Warning</b> : Do not try to access the loading tower when the transport system is accessing the loading tower . You may cause a plate transport error or crash situation if you ignore the plate loading tower alerts.
	<b>Warning</b> : Never attempt to reach the washer area while the NEO Iris is operating. You may disrupt the instrument or injure yourself. The instrument switches off power to the motors if resistance to movement is encountered.
	<b>Warning</b> : Do not remove a rack when the indicator LED is flashing red, as this can damage the pipetting system and invalidates all test results on samples in the rack. Also, do not load a rack when the LED is solid red. The barcodes are not read and the reagents or samples will not be used.
	<b>Warning</b> : Open field access is required to provide continuous access for sample/reagent loading during full system operation. Never try to access tubes or vials while their rack is still in the loading bay. Always pull their rack completely out before accessing individual positions.
	Unauthorized access to the loading bay is strictly prohibited and could injure you.
	<b>Warning</b> : It is important to keep your hands away from the pipetting area to avoid potential injury due to the moving instrument.
	<u>Warning</u> : You must not bump, knock, rest up against or otherwise come into physical contact with the centrifuge module because you may cause balancing errors or centrifuge loading and unloading issues. By causing these errors or issues, you could therefore interfere with assay processing. No items should be stored on top of the centrifuge cover including, but not limited to, sample or reagent racks, books, documents or other laboratory consumables and supplies.
	<b>Warning</b> : The liquid waste is potentially biohazardous material. Always wear protective gloves and clothing when handling the liquid waste. If any liquid waste is spilled, clean it up

lcon	Description
	immediately following the standard practice of the laboratory.
	<b>Warning</b> : Used plates and consumed liquid reagent vials contain potentially biohazardous material. Always wear protective gloves and clothing when handling used plates. If any liquid from a plate is spilled, clean it up immediately following the standard practice of the laboratory.

# **Chapter 3 – System Software Navigation**

lcon	Description
	<u>Warning</u> : You must be careful to make sure that the samples in one given rack are of a single prefix type. It is your responsibility to check this data. If you have mixed samples of different prefixes, all samples will have the same prefix added to that which you entered. Thereby, identification errors could occur if the operator fails in this tube data check.
	<u>Warning</u> : If an incorrect plate expiration date is entered, which is beyond the actual expiration date of the plate, and an assay using this plate lot is run at a date after the actual expiration date, but before the erroneously entered software expiration date, then the validity of the results is compromised, even though the results will appear legitimate. It is the operator's responsibility to enter the correct plate expiration date into the software.

# **Chapter 6 – Instrument Testing Operation**

lcon	Description
	<u>Warning</u> : Inspect all reagents and controls for the presence of foam before placing on the instrument. Do not vigorously agitate blood grouping anti-sera or controls. Shaking will produce foam in the vial that can cause the Liquid Level Detection (LLD) feature of the pipetting system to aspirate foam and/or air rather than reagent. This will produce incorrect results or an error.
	<b>Warning</b> : Before placing reagents on the instrument, you must remove the bottle caps. You are advised to remove and discard the dropper by pulling the dropper from the bulb. When you remove the reagents from the instrument for storage, you must place the caps back on the bottles. To avoid cross contamination of reagents, it is important that you place the caps on the correct bottles. Mixing caps can result in erroneous test results. Immucor-approved disposable vial caps can alternatively be used instead of replacing the saved original vial caps. Contact Technical Support for information regarding these caps.
	<b>Warning</b> : If you do not add the stirballs to the cell suspensions, the results may be invalid or incorrect. Do not touch the stirballs. You should add them directly to the cellular reagent vials using the dispenser provided. Contamination and neutralization of cellular reagents can occur if the stirballs are touched.
	stirball per vial.
	<b>Warning</b> : If you are using two or more instruments, then the specific reagent vials for each instrument must be dedicated for use on that single instrument to ensure correct reagent volume tracking. If the actual reagent volume (less than the software numeric volume) is not sufficient for the number of tests scheduled, the instrument will produce invalid results and samples will need to be rescheduled for testing.
	<b>Warning</b> : Incorrectly placing a plate in a transport frame may cause damage to the pipetting system and other modules on the instrument in addition to wasted resources.
	<b>Warning</b> : Loading a plate with an incorrect strip orientation results in invalid results and can create a biohazardous spill on the instrument. Incorrect orientation includes strips inserted upside-down in a white plate frame.

lcon	Description
	<b>Warning</b> : Loading incorrect strips in a white plate frame can produce incorrect sample results without warning. For example, loading Capture-R Ready-Screen strips in a barcoded Capture-R Select white frame is incorrect.
	<u>Warning</u> : When using Capture-R Ready-Screen or Ready-ID strips, load only the number of strips that are required for testing, as indicated in the <i>Resource Overview</i> window. Do not use Capture-R Ready-Screen or Ready-ID strips for testing that have previously traveled though the NEO Iris system on a plate frame, but were not used for testing.
	<b>Warning</b> : You must not place new plates into the tower while the transport system is accessing that module. Wait a few seconds after the transport has left the module, and then insert the new plates.
	<u>Warning</u> : You can load and store Capture strips (excluding Capture-R Ready-Screen and Ready-ID) on the instrument, outside of the storage pouch, according to the time limitations printed in the relevant package inserts. However, it is the operator's responsibility to keep track of the on-board Capture strip storage time. The instrument will provide no alerts if the on-board storage time is exceeded and will not flag results generated from such strips.
	<b>Warning</b> : Removing racks while the probes are accessing the tubes or vials from those rack results in damage to the probes and invalidated results. Wait until the LED for the rack is blinking green before removing.
	Warning: The operator must first open the loading bay dialog in order to remove sample and reagent racks.
	A continuous orange LED indicates that the samples are not required for the current processing. However, continuous orange can also indicate that outstanding tests remain for at least one sample in the rack that is not scheduled as part of the current processing.
	Warning: Prematurely removing a sample or reagent rack will invalidate results.
	<b>Warning</b> : Used plates contain potentially biohazardous material. Wear protective gloves and clothing at all times when handling used plates. If any liquid is spilled, clean it up immediately following standard laboratory practice.

# **Chapter 7 – Test Results**

lcon	Description
	Warning: Refer to Attachment 1 for Neo Operator Manual for a list of assays with
	equivocal reactions that cannot be edited.

# Chapter 10 – Maintaining the NEO Iris

lcon	Description
	<u>Warning</u> : If a maintenance action requires the manual entry of a plate ID, do not input an ID with 2 digits of alpha characters followed by numeric characters. Because the software may recognize this ID scheme as that of a sample testing plate, entry of a plate expiration date will be required. Entry of a plate expiration date is not part of the instructions for maintenance actions.
	<b>Warning</b> : Results can be adversely affected if the system liquid container is filled with anything other than PBS (e.g. de-ionized water).
	<b>Warning</b> : Filling the container with the incorrect fluid adversely affects Solid Phase Red Cell Adherence reactions.
	<u>Warning</u> : Blood samples, liquid waste, used micro-well strips, and consumed liquid reagent vials contain potentially bio-hazardous material. Always wear protective gloves and clothing when handling blood samples, liquid waste,
	used micro-well strips, or consumed liquid reagent vials. All blood samples, liquid waste, used micro-well strips, and consumed liquid reagent vials must be discarded following the standard practice of the laboratory.
	All blood products must be treated as potentially infectious. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.
	<b>Warning</b> : This maintenance task can only be performed when the NEO Iris is not active. You should only use this function for minor corrections. If a probe is badly bent, you should straighten it or replace it. Resetting a badly bent probe can result in crashes as the hardware tries to move to positions outside of the mechanical boundaries of its reach.
	<b>Warning</b> : You can open the instrument hood during this maintenance action to access the probes. When accessing the instrument deck, you must not move the pipetting turntable. Moving the turntable results in an instrument crash.
lcon	Description
------	---
	<b>Warning</b> : The main maintenance screen indicates that this action is Due if the assay successfully completes. The system does not analyze the results on the instrument, so you must determine if the volume is acceptable. You are responsible for making the corrections necessary to satisfy these requirements.
	<b>Warning</b> : Immucor requires that the quality control (QC) assay for agglutination tests is performed on the system every day, or prior to the use of agglutination reagents. The primary objective of QC is to ensure that the loaded reagents are satisfactory.
	<b>Warning</b> : Repeated QC failures in absence of a reagent related cause can be an indication of conditions that have impacted instrument system performance (i.e. pipetting, centrifugation, incubation, etc), or that an incorrect system liquid (deionized H2O rather than PBS) may have been loaded.
	Warning: For the QC3_Cell assay, the correct position for the Capture-R® Ready-Screen® (3) strip in a Capture-R® Ready-Screen® (3) white plate frame is at strip position one (1). The instrument will not check for the correct placement of the Capture-R® Ready-Screen® (3) strip in the white plate frame. You must make sure of the correct placement of this strip yourself. Incorrect placement of this single strip can result in test material being pipetted onto the pipetting station surface, an overflow error in the washer, and a failed quality control. Unused Capture-R® Ready-Screen® (3) strips should not be present in any of the other positions of the plate. It is not required that strip positions two (2) to twelve (12) of the white plate frame be occupied by any strips at all. Those positions can be left empty, but all strip positions should be activated in the software. The Strip Selection tab in the Plate Loading Tower dialog must indicate that all twelve (12) strips are available even though only the strip at position (1) is present.
	<b>Warning</b> : You should not initiate this maintenance function if an assay is scheduled to use the reader. Leaving a microplate in the reader during this action affects the analysis of results.
	<b>Warning</b> : Failure to archive results according to the configured time schedule interval, at a minimum, can result in the accumulation of excessive result data on the instrument computer which can lead to the slowing down of the speed of the computer processing and activity.

lcon	Description
	<u>Warning</u> : You must not enter the Archive tab of the Utilities window while the NEO Iris is processing, otherwise there is a risk that conversion of results could be interrupted and that results could need to be re-sent.
	<b>Warning</b> : Before you can physically empty the 20 liter system liquid container, you must remove it from the washer software. If you do not adjust the Wash Buffer dialog to reflect the buffer removal and a test is in process, the software still attempts to use the wash buffer, even when none is present. If this occurs, the associated NEO Iris tubing and pumps can run dry.
	<u>Warning</u> : The main maintenance screen indicates that this action is Due if the assay successfully completes. The system does not analyze the results on the instrument, so you must determine if the volume is acceptable. You are responsible for making the corrections necessary to satisfy these requirements.
	<u>Warning</u> : It is critical the syringe is held in perfect vertical orientation when twisting into the 3-way valve. Do not force the syringe. If there is resistance back it off and try again.
	<b>Warning</b> : Verify that you are replacing the syringe with another syringe of the same size. Replacing the syringe with a wrong-sized syringe causes pipettor inaccuracy and incorrect test results.

# Chapter 11 – Troubleshooting the NEO Iris

The following warnings are contained in this chapter:

lcon	Description
	<u>Warning</u> : During error recovery, the process control allows you to open the safety hood to access components manually. You must be cautious during manual instrument deck interventions, because while you resolve an error in one system module, other operations may remain active.
	<u>Warning</u> : After the current run is finished, you MUST shut down the software again, turn off the instrument, and then restart. If the attempt to complete the run fails, you must also shut down both the software and the instrument, and then restart. Normal operation resumes only after this complete shutdown of the instrument and the PC.

lcon	Description
	<u>Warning</u> : You should always select the <b>Abort plate</b> button (not the <b>Abort run</b> button).
	Failure to select the <b>Abort plate</b> button allows processing to continue for a plate that may
	not be valid. If you select the <b>Continue work</b> button without aborting the plate, the
	system will invalidate all results generated from the plate that suffered the pipettor error.

# Appendix A – Preparing the NEO Iris for First Use

The following warnings are contained in this appendix:

lcon	Description
	<b>Warning</b> : The power input module is being used as a separator and therefore the instrument has to be setup in a way that the appliance inlet has to be accessible from the outside in an emergency (e.g. fire) to be able to pull the power cord or plug.
4	<u>Warning</u> : Internal voltage: Always turn off the power switch and unplug the power supply before cleaning the outer surface of the instrument.
	<b>Warning</b> : <b>Potential biohazards</b> : Adequate safety precautions should be taken, as outlined in the assay's package insert. Always wear safety glasses and appropriate protective equipment, such as chemically resistant rubber gloves and apron. Dispose of waste in an approved manner. Decontaminate the instrument per the guidelines given in Chapter 10 – Maintaining the NEO Iris.
	<b>Warning</b> : <b>Unspecified use</b> : Failure to operate this equipment according to the guidelines and safeguards specified in this manual could result in a hazardous condition.
	<u>Warning</u> : Qualified personnel: Some maintenance procedures must be performed with the instrument's front cover hood elevated. Only qualified personnel—trained in the hazards involved when operating the instrument while open—are allowed to conduct these procedures, as moving parts may pose pinch or crush hazards.
	Warning: Excessive humidity: Operate the instrument on a flat surface away from excessive humidity.
	Warning: Excessive ambient light: Avoid instrument exposure to direct sunlight as it can reduce the performance range of the instrument.
	Warning: on/off switch accessibility. Do not block the right side of the instrument, as the on/off switch must be easily accessible at all times.

lcon	Description
NOTICE	Cleaning, disinfection or decontamination
	Observe the following aspects during cleaning, disinfection or
	decontamination because otherwise breakdowns or damage can be the
	result.
	• Disinfect or decontaminate components with recommended cleaning solution.
	Only use liquid cleaning, disinfection or decontamination solutions with a
	moistened cleaning tissue.
	Use only approved cleaning, disinfection or decontamination solutions and
	methods (Refer to Chapter 10: Maintaining the NEO Iris).
	Avoid cleaning, disinfection or decontamination solutions to come into contact
	with bearings and guides, as otherwise the greasy film may dissolve!
	• Do not use cleaning, disinfection or decontamination solutions in the proximity of
	circuit boards, light barriers and acrylic glass surfaces!
	• Do not pour or spray liquid cleaning, disinfection or decontamination solutions
	into the instrument.
	Do not autoclave containers and components for liquids or waste.
NOTICE	Handling of decontamination products
	Pay attention to managing the decontamination products, because they
	are harmful as indicated.
	Strictly follow the local and national provisions, legislation and laboratory
	regulations.
	Do not use bleach!
	• Do not use improper decontamination products. Use the recommended cleaning
	solution to decontaminate the instrument (Refer to Chapter 10: Maintaining the
	NEO Iris).
	• Only use decontamination solutions in accordance with the respective instructions
	for use and MSDS!
	Warning: Disposal: This instrument must be fully decontaminated prior to disposal.
	This instrument contains printed circuit boards and wiring with lead solder. Dispose of the
	instrument according to Directive 2012/19/EU, "on waste electrical and electronic
	equipment (WEEE)" or local ordinances.

# Appendix A: Preparing the NEO Iris for First Use

# In This Appendix

An Immucor representative must unpack and install the instrument because it contains sensitive electrical equipment and requires correct handling, adjustment, and testing before use. This appendix only provides the unpacking, installation, and setup information that is of interest to the operator. Contact Immucor for more information about the instrument installation.

APPENDIX A: PREPARING THE NEO IRIS FOR FIRST USE	A-1
Verifying all Parts Are Present	A-2
Environmental Conditions and General Safety Features	A-3
User Safety	A-7
Making the Connections	A-9
Software Installation	A-12
Setting Up the Instrument	A-13
Completing the Post-Installation Check	A-14
Verifying the Installation	A-15
Removal of the Instrument	A-16

# Verifying all Parts Are Present

The following is a list of all the major parts that are delivered. The Immucor representative verifies that all parts are present during the installation process.

Item	Quantity				
NEO Iris, including modules and main hood.	1				
Personal Computer (PC):					
Application software installed by manufacturer, including all applicable licenses.	1				
<b>Note:</b> Each PC belongs to a specific instrument because instrument-					
specific parameters are installed on each hard disk.					
Touch screen monitor:					
• 17" size (minimum)	1				
TFT active matrix panel					
Trackball	1				
Keyboard	1				
Barcode hand scanner	1				
Cabinet	1				
System liquid containers	2				
Common waste container	1				
Waste shuttle container	1				
Reagent ( <b>R</b> ) racks (12 position) 10ml vial	Number				
Reagent ( <b>S</b> ) racks (9 position) 10ml vial	determined by				
Reagent (T) racks (5 position) 57ml (43mm diameter) vial	specific on-site				
	requirements.				
Reagent (Z) racks (6 position) 50ml trough	1				

#### Environmental Conditions and General Safety FeaturesAppendix A: Preparing the NEO Iris for First Use

Item	Quantity				
Sample racks (16 position):					
• Tube: A Sample Rack 16–17 x 100mm					
• Tube: <b>B</b> Sample Rack 12–13 x 75–100mm					
• Tube: <b>C</b> Sample Rack	Number				
<ul> <li>Adaptor: 12–16mm outer diameter</li> </ul>	determined by				
<ul> <li>Tube: Minimum 9mm inner diameter</li> </ul>	specific on-site				
<ul> <li>Tube Height: 75mm–100mm</li> </ul>	requirements.				
<b>Note:</b> Contact Technical Support for information regarding the use of approved pediatric tubes.					
Donor: <b>D</b> Donor Rack 12mm diameter tube					
Plate carriers	15				
Dower cords and fuses	5 fuses				
	2 EU power cords				
Accessory kit	1				

# Environmental Conditions and General Safety Features

The instrument must be located in an appropriate environment. This section describes the necessary environment for the instrument and also the general safety features.

- For indoor use only.
- Installation dimensions: 86.6 inches (220 cm) width X 35.4 inches (90 cm) depth X 73.2 inches (186 cm) height. There should be an additional 22.4 inches (57 cm) of empty space in front of the cabinet to allow the cabinet doors to fully swing open. There should be an additional 6 inches (15.2 cm) of empty space above the instrument to allow raising of the hinged hood.
- Weight of the instrument: 924 lbs (420 kg). The supporting floor must be capable of bearing the load of the instrument.
- Ventilation (non-drafting environment): Installation Category per IEC 61010-1
- Pollution Degree 2: For use in areas where there is normally nonconductive pollution; however, temporary conductivity caused by condensation can occur.
- Installation category II: The instrument should be run only connected to a standard low voltage system being able to cope with potential excess voltage peaks.

• You must place the instrument on the cabinet provided. The cabinet provides ample space for storing the system liquids containers, liquid waste, PC and UPS. The cabinet also has wheels allowing for easy access to the instrument in the event that instrument service is required and it must be temporarily moved for that service work.



**Limitation:** Adequate space must be provided for ventilation and also accessing the instrument, both during operation and in the event of servicing. The cabinet has two (2) spacing pins which are 4.8 inches (12.1 cm) in length. The instrument must be positioned such that an overall clearance distance of 30 inches (76.2 cm) must be maintained between the back of the instrument and the wall. This 30 inches (76.2 cm) includes the length of the spacing pins.



**Warning:** The power input module is being used as a separator and therefore the instrument has to be setup in a way that the appliance inlet has to be accessible from the outside in an emergency (e.g. fire) to be able to pull the power cord or plug.



# A DANGER Improper connection of mains supply

Improper connection of the instrument and the peripheral devices to the mains supply can cause serious personal injury with potentially deadly consequences and material damage (e.g. fire).

- Only use grounded connection and extension cables with sufficient capacity (voltage and current) to connect the instrument and any peripheral devices to the mains power supply.
- Never remove ground connections.
- Grounding of the instrument and its peripheral devices to the same protective earth potential shall be ensured.
- The use of a multi-outlet power strip is not allowed!
- Only use power cables that fulfill the minimum requirements for this instrument.



**Warning: Internal voltage:** Always turn off the power switch and unplug the power supply before cleaning the outer surface of the instrument.



**Warning:** Potential biohazards: Adequate safety precautions should be taken, as outlined in the assay's package insert. Always wear safety glasses and appropriate protective equipment, such as chemically resistant rubber gloves and apron. Dispose of waste in an approved manner. Decontaminate the instrument per the guidelines given in **Chapter 10** – Maintaining the NEO Iris.



**Warning: Unspecified use:** Failure to operate this equipment according to the guidelines and safeguards specified in this manual could result in a hazardous condition.



**Warning:** Qualified personnel: Some maintenance procedures must be performed with the instrument's front cover hood elevated. Only qualified personnel—trained in the hazards involved when operating the instrument while open—are allowed to conduct these procedures, as moving parts may pose pinch or crush hazards.



**Limitation:** Service: The instrument should be serviced by Immucor authorized service personnel. Only qualified technical personnel should perform troubleshooting and service procedures on internal components.



**Warning: Excessive humidity:** Operate the instrument on a flat surface away from excessive humidity.



**Warning: Excessive ambient light:** Avoid instrument exposure to direct sunlight as it can reduce the performance range of the instrument.



**Limitation:** Environmental conditions: Do not expose the instrument to temperature extremes. For proper operation, ambient temperature range should remain between 18°C and 30°C. Performance may be adversely affected if temperatures fluctuate above or below this range.

- Shipping and storage temperature range: 0°C to 40°C.
- Relative humidity range: 20% to 80% (non-condensing).
- For use at altitudes up to 2,000 meters above mean sea level.



**Warning:** on/off switch accessibility. Do not block the right side of the instrument, as the on/off switch must be easily accessible at all times.

#### NOTICE

#### Cleaning, disinfection or decontamination

Observe the following aspects during cleaning, disinfection or decontamination because otherwise breakdowns or damage can be the result.

- Disinfect or decontaminate components with recommended cleaning solution.
- Only use liquid cleaning, disinfection or decontamination solutions with a moistened cleaning tissue.
- Use only approved cleaning, disinfection or decontamination solutions and methods (Refer to **Chapter 10: Maintaining the NEO Iris**).
- Avoid cleaning, disinfection or decontamination solutions to come into contact with bearings and guides, as otherwise the greasy film may dissolve!
- Do not use cleaning, disinfection or decontamination solutions in the proximity of circuit boards, light barriers and acrylic glass surfaces!
- Do not pour or spray liquid cleaning, disinfection or decontamination solutions

into the instrument.

• Do not autoclave containers and components for liquids or waste.

NOTICE Handling of decontamination products

Pay attention to managing the decontamination products, because they are harmful as indicated.

- Strictly follow the local and national provisions, legislation and laboratory regulations.
- Do not use bleach!
- Do not use improper decontamination products. Use the recommended cleaning solution to decontaminate the instrument (Refer to Chapter 10: Maintaining the NEO Iris).
- Only use decontamination solutions in accordance with the respective instructions for use and MSDS!



**Limitation:** Before using any chemicals (e.g. decontamintion solutions) except those recommended by the manufacturer, users should check with the manufacturer that chemicals will not damage the equipment or create potential hazards.



**Limitation:** Warranty: Failure to follow preventive maintenance protocols may void the warranty. Refer to Chapter 10 – Maintaining the NEO Iris.



**<u>Limitation</u>**: Changes or modifications to this unit not expressly approved by the manufacturer could **void** the user's authority to operate the equipment.

# User Safety

This device has been type tested by an independent laboratory, and it was found to meet the requirements of the following:

### **North America**

#### Underwriters Laboratories UL 61010-1:2012

"Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements"

### Canadian Standards Association CAN/CSA C22.2 Nr. 61010-1-12 (R2022)

"Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use, Part 1: General Requirements"

#### Also see other standards listed under CE Mark

### International

### IEC 61010-1:2010/AMD1:2016/COR1:2019

"Safety requirement for electrical equipment for measurement, control, and laboratory use. Part 1, General requirements"

#### Also see other standards listed under CE Mark

# **CE** Based on the testing described below and information contained herein, this instrument bears the CE mark.

#### Regulation (EU) 2017/746 In Vitro Diagnostics

Compliance with the requirements of the Directive, including Risk Assessment, Quality System Certification, and Product registration with competent authorities.

#### EN 61010-1:2010+A1:2019

"Safety requirement for electrical equipment for measurement, control and laboratory use. Part 1, General requirements"

#### EN 61010-2-020:2017

"Particular requirements for laboratory centrifuges"

#### EN 61010-2-081:2015

"Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes"

#### EN 61010-2-101:2017

"Particular requirements for in vitro diagnostic (IVD) medical equipment"

### EMC EC Directive 2014/30/CE Electromagnetic Compatibility

#### EN 61326-1 (2013)

Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1: General requirements

#### EN 61326-2-6 (2013)

Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 2 – 6: Particular requirements – In vitro diagnostic (IVD) medical equipment

(IEC 61326-2-6:2020)

#### **Emissions – CLASS A**

The system has been type tested by an independent, accredited testing laboratory and found to meet the requirements of EN 61326-2-6:2013 for Radiated Emissions and Line Conducted Emissions.

#### Immunity

The system has been type tested by an independent, accredited testing laboratory and found to meet the requirements of EN 61326-2-6:2013 for Immunity.

#### Directive 2012/19/EU Waste Electrical and Electronic Equipment

#### **Disposal Notice:**

This instrument contains printed circuit boards and wiring with lead solder. Dispose of the instrument according to Directive 2012/19/EU, "on waste electrical and electronic equipment (WEEE)."

# Making the Connections

This section outlines the connections required for the operation of the instrument, including:

- Instrument connections
- PC connections

### **Instrument Connections**

This section describes the instrument connections.

#### Descriptions

The NEO Iris connects to the main power supply through the power cord connected to the power socket on the right side of the instrument. A second power cord is connected to the back of the centrifuge module.

A Universal Serial Bus (USB) connection transfers data between the PC and the instrument. For more information about connecting the USB connection, refer to PC Connections in this appendix.

Universal Serial Bus (USB) connections connect the PC directly to the CMOS cameras in the instrument plate reader module.

The connections for the system liquid link the system liquid container in the cabinet to the washer module of the instrument. The system liquid connections are comprised of tubing for system liquid and tubing for liquid waste. A cable electronically connects the washer sub-modules, the instrument and system liquid containers.

The connections for the pipettor liquid system link the system liquid container and common waste container in the cabinet to the pipetting system of the instrument.

For correct alignment and safe placement, the rear feet of the instrument must fit into the welded metal rings on top of the cabinet.

### **PC Connections**

This section describes the PC connections.

#### Descriptions

A standard PC power cord connects the PC to the voltage supply. A USB connection transfers data between the PC and the instrument.

Connection to an LIS can be achieved by a serial connection or the network (LAN) adaptor.

The monitor, trackball and keyboard (via the barcode hand scanner) are connected in the same manner as a standard PC. The touch screen electronics of the monitor are connected to the USB connection.

A standard PC power cord connects the monitor to the voltage supply.

The power cords of the monitor, PC and instrument are connected to the UPS.

#### Using an Uninterruptible Power Supply (UPS)

The NEO Iris is a computer-controlled instrument that electronically stores both test requests and test results. You can lose the current run data (but not the data that has previously been generated) if the power supply to the system is interrupted. Using an Uninterruptible Power Supply (UPS) with an incorporated power conditioner is required to operate the instrument. Approximately 15 minutes of battery power is available from the UPS to run the instrument in the absence of wall outlet power.

# Software Installation

The software for the instrument is preinstalled on the PC.



**Limitation**: Do not install additional software on the PC. This will void your warranty and service contract. Third-party software can interfere with the controlling software and result in a loss of sample data.

# Setting Up the Instrument

The instrument is dedicated for use with Immucor reagents. Nevertheless, you can configure several functions according to the specific needs of the laboratory. You can configure the following parameters:

- Operator name and access rights
- Operator password
- Frequency of archiving results

For more information about configuration, refer to **Chapter 4 – Security**.

# Completing the Post-Installation Check

After the Immucor representative has finished the installation, you must check all of the components of the instrument for proper functioning. To check the components, you must follow the procedures described in **Chapter 10 – Maintaining the NEO Iris**.

# Verifying the Installation

An installation and operational qualification process is completed by Immucor after the instrument has been installed.

The installation and operational qualification process confirms that:

- The operating environmental conditions are met.
- The Immucor documentation is supplied.
- The mechanical safety features are active.
- All calibration results are acceptable.
- The necessary instrument modules function correctly.
- The assays that have been configured for laboratory requirements are correctly prepared and interpreted.

# Removal of the Instrument

This section describes the steps to be taken prior to:

- Shipment of the instrument back to the manufacturer.
- Waste disposal of the instrument.

#### Shipment to Manufacturer

If it is necessary to remove the instrument from the laboratory and ship it back to the manufacturer, an Immucor representative must repackage it into the original packaging.



**Limitation: Warranty:** You must save all of the original packaging materials. If you need to ship the instrument to Immucor for repair or replacement, the original packing must be used. Other forms of commercially available packing are not recommended and can void your warranty.



**Limitation**: You must decontaminate the instrument before it is repackaged. Follow the procedures in **Chapter 10 – Maintaining the NEO Iris** to decontaminate the instrument. The responsible laboratory authorities must state on the shipping documents that the instrument has been decontaminated and provide a signature for verification. Immucor reserves the right to refuse to accept an instrument that does not have such a decontamination statement.

#### **Disposal to Waste Facilities**

If it is necessary to dispose of the instrument to waste facilities, you must decontaminate the instrument before it is discarded. Follow the procedures in **Chapter 10 – Maintaining the NEO Iris** to decontaminate the instrument.



Warning: Disposal: This instrument must be fully decontaminated prior to disposal.

This instrument contains printed circuit boards and wiring with lead solder. Dispose of the instrument according to Directive 2012/19/EU, "on waste electrical and electronic equipment (WEEE)" or local ordinances.

# **Appendix B: Maintenance Records**

# In This Appendix

This appendix provides master forms of the NEO Iris maintenance records.

Use copies of these master forms to record results of the maintenance procedures as described in **Chapter 10 – Maintaining the NEO Iris**.

APPENDIX B: MAINTENANCE RECORDS	B-1
Maintenance Forms	В-2
NEO Iris Daily Maintenance Record	B-3
NEO Iris Weekly Maintenance Record	В-4
NEO Monthly Maintenance Record	B-5
NEO Iris PipTest Chart: Precision and Accuracy	В-6
NEO Iris Residual Volume Calculation Record	B-7

# Maintenance Forms

The following forms are in this section:

- NEO Iris Daily Maintenance Record
- NEO Iris Weekly Maintenance Record
- NEO Iris Monthly Maintenance Record
- NEO Iris PipTest Chart: Precision and Accuracy
- NEO Iris Residual Volume Calculation Record

# **NEO Iris Daily Maintenance Record**

**Instructions for use:** Operator(s) must perform the maintenance requirement at the specified time. Ink initials of the operator in the designated box(es) signify successful completion of the task(s). Refer to **Chapter 10 – Maintaining the NEO Iris** for written information about daily maintenance requirements.

	FACILITY								MONTH/YEAR											INSTRUMENT #												
MAINTENANCE REQUIREM	ENT																															
	DATE	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	8 19	20	21	22	23	24	25	26	27	28	29	30	31
Hardware and Fluids:																																
Perform a complete NEO Iris mode shutdown followed immediately by powering both items on.	ule and PC /																															
Check and refill system liquid cont necessary.	ainer if																															
Empty waste container.																																
Clean and initialize: check that liqu coming outside the probes and th overflow occurs.	iid is at no																															
Check the pipettor reference.																																
Pipettor Self Check																																
0	perator																															

# **NEO Iris Weekly Maintenance Record**

**Instructions for use:** Operator(s) must perform the maintenance requirement at the specified time. The date of each weekly activity must be entered in the table. Ink initials of the operator in the designated box(es) signify successful completion of the task(s). Refer to **Chapter 10 – Maintaining the NEO Iris** for written information about weekly maintenance requirements.

FACILITY		INSTRUMENT #								
MONTH/YEAR:										
MAINTENANCE	REQUIREMENT									
		Operator	Date							
Perform the act	ions below in the numbered sequence:									
1. Power <u>off</u> th	ne NEO Iris module ( <b>only</b> the module <u>and not</u> the PC).									
2. Fill probe rin	se cup and waste cup of both pipettor wash towers with recomm	nended cleaning solution;								
and expose f	for prescribed contact time.									
3. Common wa	aste Container:									
(a) Empty and ri	nse with the recommended cleaning solution. (b) Clean inside of	lid with alcohol wipes or								
recommended cl	leaning solution. (c) Empty and rinse with deionized water. (d) Er	npty and reconnect the								
container.										
4. Inspect syring	ges for leakage and check tightness of barrels.									
5. Power on th	e NEO Iris module and allow initialization to occur.									
N N	ote: the instrument must be powered back on for a minimum of	thirty (30) minutes to allow								
th	the reader lamp to warm up before using the reader module to read plates.									
6. Reader: Clea	6. Reader: Clean the mirror and Take new Flat-field images									
7. Perform resu	Ilts Archive.									

### **NEO Monthly Maintenance Record**

**Instructions for use:** Operator(s) must perform the maintenance requirement at the specified time. The date of each monthly activity must be entered in the table. Ink initials of the operator in the designated box(es) signify successful completion of the task(s). Refer to **Chapter 10 – Maintaining the NEO Iris** for detailed written information about monthly maintenance requirements.

FACILITY		YEAR	MONTH		INSTRUMENT #	
MAINTENAN	CE REQUIREMENT					
Operator Notation: Operator Date						Date
Software: Cle	ar Test Data.					
Run the <b>DECC</b> following spe	<b>DNTAMINATION</b> procedure ( <b>Decontaminate tu</b> cific steps:					
After he <b>DB_CK</b> assay is finished, verify that liquid has collected into the collection container from the blue tubing.						
After the <b>Prin</b>	<b>ne3</b> assay is finished, verify that there is no visib	late.				

# NEO Iris PipTest Chart: Precision and Accuracy

Refer to Chapter 10 - Maintaining the NEO Iris for written information about the Pipettor Verification Test (PipTest)

Strip	1	2	3	4	5	6	7	8	9	10	11	12
Beginning Weight												
Ending Weight												
Actual Volume												
Expected Volume	0.18-0.22	0.36-0.44	0.54-0.66	0.72-0.88	0.36-0.44	0.36-0.44	0.36-0.44	0.36-0.44	0.75-0.92	0.75-0.92	0.75-0.92	0.75-0.92
Outcomes acceptable ? Check one (1) box only (Yes or No).	Yes											
	No											
Instrument Number:												
Name of Site:												
Operator:												
Date:												

### **NEO Iris Residual Volume Calculation Record**

Instructions for use: Operator(s) must perform this maintenance action on an as needed basis only at the request of Technical Support.

Ink initials of the operator in the designated box signify completion of the task. Refer to **Chapter 10 – Maintaining the NEO Iris** for written information about the residual volume action.

FACILITY	INSTRUMENT serial number	
----------	--------------------------	--

NEO Iris Residual Volume	Item descriptor				
Operator notation:	Operator	Date	Weight	Scale Serial Number	
<b>A.</b> Record operator initials and date of performance:					
<b>B.</b> Record weight of plate <u>before</u> test process:			g		
<b>C.</b> Record weight of plate <u>after</u> test process:			g		
D. Subtract weight <u>before</u> test from weight <u>after</u> test:			g		
E. Record serial number of electronic scale:					

# **Appendix C: Hardware Technical Data**

### In This Appendix

This appendix contains a summary of both general technical data about the instrument and module specific data.

APPENDIX C:	HARDWARE TECHNICAL DATA	:-1
Hardware Te	echnical Data	:-2

# Hardware Technical Data

# **General Data**

This section describes the general technical data for the instrument.

Area	Requirement	
Power Requirements		
Voltage	100 – 240 V	
Frequency	50 – 60 Hz	
Current consumption	4 – 1.7 A	
Mean power consumption	250 W	
Fuse	250 VAC, T4 AH	
Dimensions		
Instrument with cabinet	<ul> <li>Width (without monitor track): 63 inches (160 cm)</li> <li>Width (with monitor track): 86.6 inches (220 cm)</li> <li>Depth: 35.4 inches (90 cm)</li> <li>Height: 73.2 inches (186 cm)</li> </ul>	
Weight		
Instrument	926 lbs (420 kg)	
Connections		
LIS to PC	Serial RS232 or LAN (network)	
PC to instrument	USB	
Internal instrument	CAN bus	

# Personal Computer (PC)

This section describes the hardware technical data for the instrument PC and its software.

Area	Requirement
Hardware	
Processor	Intel Core 2 Duo
Memory	Greater than or equal to 1 GB
Drives	1 DVD drive
Ports	USB (2 front and 6 rear)
	1 parallel
Touch-screen monitor	TFT active matrix panel
Keyboard	QWERTY, AZERTY, or QWERZ (depending upon local requirements)
Handheld barcode scanner	Keyboard wedge type
Trackball	PS/2 Microsoft compatible
Software	
Operating System	Windows 7
Software version	Greater than or equal to 1.6.11
Data storage capacity	Greater than or equal to 80 GB, capable of storing approximately 225,000 results.

### Cabinet

This section describes the hardware technical data for the instrument cabinet.

Area	Description
Doors	Three Doors
System liquid	
Common waste container	1
Number of system liquid containers	2

### Plate Loading Tower

This section describes the hardware technical data for the instrument plate loading tower.

Area	Description
	• Width: 7.2 inches (18.4 cm)
Dimensions	• Depth: 4.4 inches (11.2 cm)
	• Height: 15 inches (38.2 cm)
Positions	15

### Transport System

This section describes the hardware technical data for the instrument transport system, including the plate barcode scanner.

Area	Description
	• Width: 45.3 inches (115 cm)
Dimensions	• Depth: 5.3 inches (13.5 cm)
	• Height: 16.9 inches (43 cm)
Transportation Time	Maximum 1 minute

Plate Barcode Scanner	
	• Width: 3.1inches (8 cm)
Dimensions	• Depth: 2 inches (5 cm)
	• Height: 0.9 inches (2.25 cm)
Connection	RS232
Barcode type	Code 128
Optimum reading distance	2 inches (5 cm)

# 14-lane and 5-lane Bays

This section describes the hardware technical data for the instrument 14-lane and 5-lane bays.

Area	Description		
14-lane bay			
	• Width: 14.8 inches (37.6 cm)		
Dimensions	• Depth: 15.9 inches (40.3 cm)		
	• Height: 4.8 inches (12.2 cm)		
Lanes	14		
5-lane bay			
	• Width: 6.4 inches (16.2 cm)		
Dimensions	• Depth: 12.8 inches (32.4 cm)		
	• Height: 4.8 inches (12.2 cm)		
Lanes	5		
Sample/reagent vial barcode scanner			
	• Codabar		
	• Code 128		
Barcode types	• Code 39		
	Interleaved 2 of 5		
	ISBT 128 (concatenated barcodes are not supported)		
	Module width must be above 0.167 mm for the narrowest		
	barcode feature.		
Barcode quality	Ratio: greater than or equal to 1:2.		
	Print quality must be Class A or B ANSI X3.182 or DIN EN		
	1635.		
Maximum fraction of reading failures	1:104		
Maximum fraction of reading error	1 : 10 <sup>7</sup>		

# Pipetting System

This section describes the hardware technical data for the instrument pipetting system.

Area	Description	
Pipettors		
Wash station dimensions	Width: 1.3 inches (3.4 cm)	
	Depth: 2.4 inches (6 cm)	
	Height: 5.9 inches (14.9 cm)	
Modes of function	Aspirating	
	Dispensing	
	Multi dispensing	
	Priming	
	Washing (steel probes)	
Channels	1x4 probes, 1x1 probes	
Volume range	Stainless steel probes: 1–1000 μl	
	Definable in increments of 1µl	
Dispensing accuracy	<2% for 100 to 500 μl	
	<5% for 25 μl	
Aspirating accuracy	<5% for 10 μl	
Liquid viscosity	Maximum 10 mPa	
System liquid container capacity	20L system liquid container connected to 10 liter refill container	
Pump Rack		
Dimensions	Width: 14.9 inches (37.9 cm)	
	Depth: 8.7 inches (22 cm)	
	Height: 13.8 inches (35 cm)	

Area	Description
Probe washing pumps	Membrane type pumps
Syringe pumps	Sample: 500 μl
	Reagent: 1000 μl
## Incubator

This section describes the hardware technical data for the instrument incubator.

Area	Description
	Width: 8.6 inches (21.9 cm)
Dimensions	Depth: 7.3 inches (18.5 cm)
	Height: 14.9 inches (37.8 cm)
T	38°C – 40°C zone
Temperature range	Ambient zone
Setting resolution	1°C
E en em tien	< 2 μl per well
Evaporation	(Uncovered plates at 37°C for 1 hour with 100 $\mu$ l distilled water per well).
Incubation Time	
Range	5-1440 minutes
Resolution	1 min

### Washer

This section describes the hardware technical data for the instrument washer.

Area	Description
	Width: 8.7 inches (22 cm)
Dimensions	Depth: 7.5 inches (19 cm)
	Height: 14.8 inches (37.5 cm)
Modes of operation	Plate mode and strip mode
Aspiration time	Programmable
Dispense volume range	50 to 2500 $\mu l$ per well with direct overfill aspiration.
Viscosity of initial well contents	Up to 10 mPa
Maximum number of wash steps	8

# Centrifuge

This section describes the hardware technical data for the instrument centrifuge.

Power Requirements	
Voltage	100 – 240 V
Frequency	50 – 60 Hz
Current consumption	4 – 1.7 A
Power consumption	200 W
Fuse	250 VAC, T4 AH
Area	Description

	Width: 21.38. inches (54.3 cm)
Dimensions	Depth: 15.28 inches (38.8 cm)
	Height: 22.76 inches (57.8 cm)
Weight	95 lbs (43 kg)
Modes of operation	Spinning Shaking
g-force	2–1200 RCF

Connections	
Description of connector	Connected to
Centrifuge Loading Unit	Loading Unit at the left side of the NEO Iris
	instrument
RS232	PC (For service use only)
CAN Connector	NEO Iris instrument
Power Supply	Separate power connector

Instrument setup (connections) must be performed by trained service personnel only (refer to **Limitation:** Service and Warning: Qualified personnel in Chapter 12: Limitations of Use and Warnings and in Appendix A)

### Camera Reader

This section describes the hardware technical data for the instrument camera reader.

Area	Description
	Width: 8.7 inches (22 cm)
Dimensions	Depth: 7.5 inches (19 cm)
	Height: 14.8 inches (37.5 cm)
Modes of operation	Image reading

Resolution	1280 x 1024 Plate image
Resolution	140 x 140 Well image

# **Glossary of Terms**

### In This Glossary

This glossary contains the definitions of the terms used in the NEO Iris Operator Manual regarding the NEO Iris hardware and software.

# Glossary

The table below lists the terms found in the NEO Iris Operator Manual in alphabetical order with their corresponding definitions.

Term	Definition
14-lane and 5-lane bays	The 14-lane and 5-lane bays are areas of the NEO Iris deck equipped to contain the blood sample tubes and reagent vials required for pipetting. The 14-lane bay is designed to accommodate both sample tubes and reagent vials. The 5-lane bay is designed to only accommodate reagent vials, and not sample tubes.
Assay Profile	Assay Profile buttons are site-customized buttons that can be built to combine two or more assays. Profiles are a group of assays that can be selected together and therefore make assay selection more efficient.
Barcode	A barcode is a series of varying width vertical lines (called bars) and spaces. There are different combinations of the bars and spaces, which represent different characters. Unique barcode character sequences are used to uniquely identify an item that the barcode is affixed to.
Cabinet	The NEO Iris cabinet is used as a surface on which to place the NEO Iris, as well as a compartment in which to store external system components, such as the computer.
Centrifuge	The on-board centrifuge is located at the left-hand side of the NEO Iris and is used to spin or shake microplate strips.
Common Waste Container	The common waste container is located in the cabinet and is used as a system collection vessel for liquid biohazardous waste.
Emergency Stop Button	The Emergency Stop button is located in the middle of the on-screen Status Bar and if pressed, will stop all NEO Iris actions.
Event Log	The Event Log is located to the left of the Emergency Stop button on the Status Bar and displays a list of event messages that have been generated since the last database clean-up.
Export Protocol	The Export Protocol describes the computer file structure used in the transfer of electronic data from the computer to a Laboratory Information System (LIS).

Term	Definition
Fuse	A fuse is an electrical safety device which will blow or melt if the current flowing through it exceeds a specified value.
Incubator	The incubator provides the necessary temperature environment for the incubation steps of the assays performed on the NEO Iris platform.
Initialization	Initialization is a NEO Iris process whereby the whole of the NEO Iris system is reset. This process includes returning all modules to their home positions.
Instrument Settings	The instrument settings piece of the NEO Iris software is the software portal where the operator user name, password, and user rights are defined.
Light Emitting Diode (LED)	LEDs are diodes that contain a semi-conductor chip that emits light within a very narrow frequency range when connected in a circuit. LEDs that emit different colors are made of different semi-conductor materials and require different energies to light them.
Login	The Login button displayed on the Main Menu Bar of the NEO Iris software screen is the portal button to access the Login screen that allows an operator to log into the NEO Iris software.
Machine Monitor	The Machine Monitor is a graphical user interface displayed on the NEO Iris software screen that reflects both the deck structure of the NEO Iris and the pattern of the workflow.
Main Menu Bar	The Main Menu Bar is displayed on the top of the NEO Iris software screen and is composed of several portal buttons that, if pressed, allow operator access to various NEO Iris software functional areas.
Maintenance	Maintenance actions for the NEO Iris are performed on a scheduled basis to verify that specific modules of the NEO Iris are functioning at the required specifications.
Manifold	The manifold is the detachable wash head, which when correctly secured in the NEO Iris microplate washing module, is used to physically dispense clean wash fluid into test wells requiring a wash step and also remove contaminated post-wash fluid out of the test wells.

Term	Definition
NEO Iris	The NEO Iris is an immunohematology device designed as a platform for in vitro diagnostic blood sample testing.
Passwords and User Access Rights	The NEO Iris employs user names, passwords, and access rights to restrict access to various levels of the NEO Iris software to authorized operators only. Unauthorized operators are refused access.
Pipetting System	The pipetting system aspirates liquids from a defined source and deposits them in a defined destination.
Plate Carrier	All microplates are transported on plate carriers. The plate carrier is a plastic frame containing springs and opposed mechanical stops on the inside of the frame to hold the microplate firmly and correctly in place.
Plate List	The Plate List displays data regarding all microplates that have been processed. The operator can also use the plate list to search for, view, and delete microplate data.
Plate Loading Tower	Plates are loaded into and removed from the NEO Iris via the plate loading tower, which is located on the left hand side of the NEO Iris deck.
Pump	A pump is a machine or device for raising, compressing, or transferring fluids.
Racks	Racks are used on the NEO Iris for the loading and unloading of sample tubes and reagent vials into and out of the 14-lane and 5-lane bays on the NEO Iris deck.
Regions of Interest (ROI)	ROIs are predetermined coordinates on a camera image of a microplate well that are analyzed through an algorithm to determine the result of that well.
Report	A NEO Iris report is printed blood sample data organized in a logical sequence, such as by numeric sequence or by calendar sequence.
Schedule	A schedule is a plan for performing work or achieving an objective, specifying the order and allotted time for each part. The NEO Iris schedule is the sequence of events that the NEO Iris must perform to complete the requested assays.
Shutdown	Shutdown is a cessation of operations or activity. The NEO Iris shutdown procedure guides the operator through the necessary steps to close down the NEO Iris.

Term	Definition
Start Run Assistant	The Start Run Assistant button, located on the Main Menu Bar, is used to access the Start Run Assistant. The Start Run Assistant is an intuitive guide for the operator to start assay runs.
Trackball	A trackball is a computer screen pointing device that does not require movement on a pad, but uses manual rotation of a ball around a fixed coordinate to move the visible pointer across the computer screen.
Transport System	The transport system moves microplates (each plate being held in a plate carrier) between the different NEO Iris modules located across the NEO Iris deck.
Touch-Screen Monitor	The NEO Iris monitor screen is pressure sensitive such that, by pressing the screen at a designated coordinate, the operator can initiate a designated NEO Iris action.
Upload Worklist	The Upload Worklist function is the process used to upload a worklist from an LIS or external computer to the NEO Iris computer.
UPS (or UPM)	The NEO Iris computer is connected to an Uninterruptible Power Supply (UPS) or Uninterruptible Power Manager (UPM) with an integral power conditioner. The UPS (or UPM) provides a consistent source of electricity (without excessive electrical noise output voltage) and also provides short periods of backup power in the event of power loss in the vicinity of the NEO Iris.
Utilities	The NEO Iris Utilities function allows the operator to view event logs and statistics as well as archive data.
Wash Stations (four probes and single probe)	Wash stations for both the four probes and single probe arms are used to clean the probes after a pipetting procedure. This cleaning prevents cross- contamination due to carry-over of blood samples and/or reagents. The wash stations are also used to prime the liquid system prior to its first use.
Washer (Microplate)	The NEO Iris microplate washer is used to physically dispense clean wash fluid into test wells requiring a wash step and also remove contaminated post-wash fluid out of the test wells. The washer uses a predetermined volume of wash fluid per well.

Term	Definition
Worklist Editor	The NEO Iris Worklist Editor function allows the operator to view, edit, and request assays for blood samples already stored in the NEO Iris software.
XYZ Movement	The transport system allows movement of the microplate and its plate carrier in the X, Y, and Z direction. The X-direction is defined as being horizontal left to right NEO Iris. The Y-direction is defined as being horizontal front to back. The Z-direction is defined as being vertical top to bottom.

# Attachment I: NEO Iris Operator Manual

### In This Attachment

This attachment describes the reagents and cutoff values used for the NEO Iris assays. It also includes basic assay procedural steps.

ATTACHMENT I: NEO IRIS OPERATOR MANUAL	I-1
Copyrights and Disclaimers	
Sample Requirements	
Assay Descriptions	
Assay Cutoffs	I-24
Assay Reagent Component Grid	I-48
Assay Procedural Steps	I-51
Test Results and Interpretation	I-82

# Copyrights and Disclaimers

©2016, Immucor, Inc. NEO Iris™ is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment 1 for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrol-de@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

### **Document Revision History**

Date	Version	Page	Description
MAY 2016	NEO Iris_EU- 001-100 (A-I)	N/A	First Version of Attachment 1 for NEO Iris Operating Manual
JUL2016	NEO Iris_EU- 001-101 (A-I)	1-4	Update of the limitations on interfering substances for Lipemia and Icterus (Levels defined)
JUN2017	NEO Iris_EU- 001-102 (A-I)	Attachment 1: NEO Iris Operator Manual	Manufacturer new address and Document Revision History List of assays updated Change the Anti-A positive cutoff range for samples from 70 to 100 to 58 to 100 for ABO assays
NOV2018	NEO Iris_EU- 001-103 (A-I)	1-12, 133, 134, 153, 157	Capture-CMV included
SEP2021	NEO Iris_EU- 002-103 (A-I)	1-29	Revised Cut-Offs for the Crossmatch assays (XMatch_E and XMatch_S).
JAN2022	NEO Iris_EU- 003-103 (A-I)	1-9	Anti-D NOVACLONE is used in the vertical ABDCHECK assay
ОСТ2022	NEO Iris_EU- 004-103 (A-I)	I-7,I-8, I-18,I-19	immuClone Rh-Hr Control is used in the horizontal ABDLONG and BABY_BG assays

# Sample Requirements



**Limitation:** At least 250 µl of packed red blood cells need to be present in a sample tube to ensure that the probe picks up red blood cells and not plasma (only for those assays that require red blood cells).

At least 500  $\mu$ l of plasma or serum needs to be present in a sample tube to ensure that the probe picks up plasma or serum.



**Limitation:** Red blood cell samples only collected with EDTA anticoagulants can be tested on the instrument. The use of other anticoagulants for testing on the instrument must be verified in the corresponding reagent Instructions for Use. Serum samples can also be tested on the instrument for tests that do not require red cells. Samples obtained from tubes containing neutral gel separators may produce falsely positive results and should therefore not be tested on the instrument.



**Limitation:** The dead volume, or the depth to which probes cannot reach, is 200  $\mu$ l of sample.



**Note:** Red blood cell donor unit segments can be tested on the instrument using the forward ABO and Rh blood grouping assays and Crossmatch. You must remove segments from blood bags and cut them so that the contents can be dispensed into a test tube with a diameter of 12 mm. It is recommended that you label the test tube with a unit ID barcoded label taken from the correct blood bag. You must then insert the test tube into the Donor rack for testing on the instrument.



**Limitation:** When loading a test tube, containing red blood cell donor unit segment contents, into a Donor rack for compatibility testing with the Crossmatch assay, the contents of the test tube must not exceed a liquid height of 3.3cm (1.3 inches), otherwise the clot detection mechanism can be falsely alarmed.



**Limitation:** Samples that exhibit excessive hemolysis should not be tested on the instrument.Samples that exhibit a hemolysis grade of 3+ or greater must not be tested on the instrument, because they may generate erroneous results. Refer to the photograph below for a hemolysis grade of 3+ (color guide). For assays using Capture-R® Select, do not use hemolyzed samples for creating a monolayer. Fragmented red blood cell membranes will interfere with monolayer formation.





**Limitation:** Samples that exhibit excessive lipemia should not be tested on the instrument. Samples that exhibit a dosage above 600 mg/dL of lipids must not be tested on the instrument because they may generate erroneous results.



**Limitation** :Samples that exhibit excessive icterus should not be tested on the instrument. Samples that exhibit a dosage above 30 mg/dL of bilirubin must not be tested on the instrument because they may generate erroneous results.



**Limitation:** Instrument studies have demonstrated that the probe washes are sufficient to prevent carry-over of samples that have an antibody titer of up to 1024 (typical of those samples encountered in a blood bank setting). It is important to note that these studies apply only to routine patient or donor samples. Sample material that is provided in surveys has been shown to occasionally cause carry-over at lower titers (a result of the manufacturing process for these materials).



**Limitation:** For assays using Capture-R® Ready-Screen® and Capture-R® Ready-ID® microplates, avoid testing previously frozen specimens that have undergone multiple freeze-thaw cycles, as these may give erroneous or inconsistent test results.

# Assay Descriptions

# List of Assays

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Horizontal/12 wells per strip	
Forward ABO	ABDCHECK	1. NOVACLONE Anti-A	Untreated microplates
Blood Grouping		2. NOVACLONE Anti-B	
		3. NOVACLONE Anti-D	
		4. NOVACLONE Diluent Control	
Forward ABO	ABDCHECK2	1. immuClone Anti-A	Untreated microplates
Blood Grouping		2. immuClone Anti-B	
		3. immuClone Anti-D rapid	
		4. immuClone Rh-Hr Control	
Forward and	ABORH	1. NOVACLONE Anti-A	Untreated microplates
reverse ABO		2. NOVACLONE Anti-B	
Blood Grouping		3. Referencells A1	
		4. Referencells B	
		5. NOVACLONE Anti-D	
		6. NOVACLONE Diluent	
		Control	
Forward and	ABORH2	1. immuClone Anti-A	Untreated microplates
reverse ABO		2. immuClone Anti-B	
Blood Grouping		3. Referencells A1	
		4. Referencells B	
		5. immuClone Anti-D rapid	
		6. immuClone Rh-Hr Control	

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Horizontal/12 wells per strip	
Forward and	ABDLONG	1. NOVACLONE Anti-A	Untreated microplates
reverse ABO		2. NOVACLONE Anti-B	
Blood Grouping		3. NOVACLONE Anti-AB	
		4. Referencells A1	
		5. Referencells A2	
		6. Referencells B	
		7. Referencells O	
		8. Autocontrol	
		9. NOVACLONE Anti-D	
		10. immuClone Anti-D rapid	
		11. immuClone Rh-Hr Control	
Forward and	ABDLONG2	1. immuClone Anti-A	Untreated microplates
reverse ABO		2. immuClone Anti-B	
Blood Grouping		3. immuClone Anti-AB	
		4. Referencells A1	
		5. Referencells A2	
		6. Referencells B	
		7. Referencells O	
		8. Autocontrol	
		9. NOVACLONE Anti-D	
		10. immuClone Anti-D rapid	
		11. immuClone Rh-Hr Control	

|--|

	Blood Grouping	Horizontal/12 wells per strip	
Forward ABO	BABY_BG	1. immuClone Anti-A	Untreated microplates
Blood Grouping		2. NOVACLONE Anti-A	
		3. immuClone Anti-B	
		4. NOVACLONE Anti-B	
		5. immuClone Anti-AB	
		6. NOVACLONE Anti-AB	
		7. NOVACLONE Anti-D	
		8. immuClone Anti-D rapid	
		9. immuClone Rh-Hr	
		Control	
Subgroup "A"	A_subg	1. Anti-A1 (Lectin)	Untreated microplates
testing		2. immuClone Anti-H	
		3. immuClone Rh-Hr	
		Control	
Forward ABO	AB_CTR	1. NOVACLONE Anti-A	Untreated microplates
Blood Grouping		2. NOVACLONE Anti-B	
		3. NOVACLONE Diluent	
		Control	
Forward ABO	AB_CTR2	1. immuClone Anti-A	Untreated microplates
Blood Grouping		2. immuClone Anti-B	
		3. immuClone Rh-Hr	
		Control	

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Horizontal/12 wells per strip	
Rh Blood Group	RHFORMEL	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-C(2)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-c(2)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-E(2)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Anti-K(1)</li> <li>Automated Anti-Kell</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
Rh Blood Group Rh Blood Group	PHENO16 AG_K	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Anti-K(1)</li> <li>immuClone Rh-Hr Control</li> <li>immuClone Anti-K(1)</li> </ol>	Untreated microplates Untreated microplates
(Specific Kell testing)		2. immuClone Rh-Hr Control	

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Horizontal/12 wells per strip	
Rh Blood Group (Specific Kell testing)	Kell	<ol> <li>immuClone Anti-K(1)</li> <li>Automated Anti-Kell</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
Rare phenotype testing	AG_M	<ol> <li>immuClone Anti-M</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
Rare phenotype testing	Cw_lgM	<ol> <li>immuClone Anti-Cw</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Vertical/8 wells per strip	
Forward ABO Blood Grouping	ABDCHECK	1. NOVACLONE Diluent Control	Untreated microplates
		2. NOVACLONE Anti-A	
		3. NOVACLONE Anti-B	
		4. NOVACLONE Anti-D	

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Vertical/8 wells per strip	
Forward ABO Blood Grouping	ABDCHECK2	<ol> <li>immuClone Rh-Hr Control</li> <li>immuClone Anti-A</li> <li>immuClone Anti-B</li> <li>immuClone Anti-D rapid</li> </ol>	Untreated microplates
Forward and reverse ABO Blood Grouping	ABDFULL	<ol> <li>NOVACLONE Diluent Control</li> <li>NOVACLONE Anti-A</li> <li>NOVACLONE Anti-B</li> <li>NOVACLONE Anti-AB</li> <li>immuClone Anti-D rapid</li> <li>NOVACLONE Anti-D</li> <li>Referencells A1</li> <li>Referencells B</li> </ol>	Untreated microplates
Forward and reverse ABO Blood Grouping	ABDFULL2	<ol> <li>immuClone Rh-Hr Control</li> <li>immuClone Anti-A</li> <li>immuClone Anti-B</li> <li>immuClone Anti-AB</li> <li>immuClone Anti-D rapid</li> <li>NOVACLONE Anti-D</li> <li>Referencells A1</li> <li>Referencells B</li> </ol>	Untreated microplates

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Vertical/8 wells per strip	
Forward and reverse ABO Blood Grouping	ABOD12_I	<ol> <li>immuClone Rh-Hr Control</li> <li>immuClone Anti-A</li> <li>immuClone Anti-B</li> <li>immuClone Anti-AB</li> <li>immuClone Anti-D rapid</li> </ol>	Untreated microplates
		<ol> <li>6. Referencells A1</li> <li>7. Referencells B</li> </ol>	
Forward and reverse ABO Blood Grouping	ABOD12NC	<ol> <li>NOVACLONE Diluent Control</li> <li>NOVACLONE Anti-A</li> <li>NOVACLONE Anti-B</li> <li>NOVACLONE Anti-AB</li> <li>NOVACLONE Anti-D</li> <li>Referencells A1</li> <li>Referencells B</li> </ol>	Untreated microplates
Forward ABO Blood Grouping	Baby_BG	<ol> <li>NOVACLONE Diluent Control</li> <li>NOVACLONE Anti-A</li> <li>NOVACLONE Anti-B</li> <li>NOVACLONE Anti-AB</li> <li>immuClone Anti-D rapid</li> <li>NOVACLONE Anti-D</li> </ol>	Untreated microplates

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Vertical/8 wells per strip	
Forward ABO Blood Grouping	Baby_BG2	<ol> <li>immuClone Rh-Hr Control</li> <li>immuClone Anti-A</li> <li>immuClone Anti-B</li> <li>immuClone Anti-AB</li> <li>immuClone Anti-D rapid</li> <li>NOVACLONE Anti-D</li> </ol>	Untreated microplates

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Vertical/8 wells per strip	
Forward ABO Blood Grouping	ABODBB_I	<ol> <li>immuClone Rh-Hr Control</li> <li>immuClone Anti-A</li> <li>immuClone Anti-B</li> <li>immuClone Anti-AB</li> <li>immuClone Anti-D rapid</li> </ol>	Untreated microplates
Subgroup "A" testing	A_subg	<ol> <li>Anti-A1 (Lectin)</li> <li>immuClone Anti-H</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
Substance "H" testing	AG_H	<ol> <li>immuClone Anti-H</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Vertical/8 wells per strip	
Rh Blood Group	PHENO12	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Anti-K(1)</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
Rh Blood Group	PHENO	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-e(1)</li> <li>Automated Anti-Kell</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Blood Grouping	Vertical/8 wells per strip	
Rh Blood Group (Specific Kell testing)	AG_K	<ol> <li>Automated Anti-Kell</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
Rare phenotype testing	AG_M	<ol> <li>immuClone Anti-M</li> <li>immuClone Rh-Hr</li> <li>Control</li> </ol>	Untreated microplates

Rare phenotype	Cw_lgM	1. immuClone Anti-Cw	Untreated microplates
testing		2. immuClone Rh-Hr	
		Control	

Assay Description	Assay Short Name	Used Reagents	Microplates used
	QC		
QC	QC	According to customer specific requirements	Untreated microplates

Assay Description	Assay Short Name	Generic name	Microplates used	
	Capture-R Ready Screen			
RBC antibody screen (pooled cells)	Pool_Cell	Capture-R RS Pooled Cells	Capture-R RS Pooled	
RBC antibody screen (2-cell)	2_Cell	Capture-R RS 2 Cells	Capture-R RS I and II	
RBC antibody screen (3-cell)	3_Cell	Capture-R RS 3 Cells	Capture-R RS 3	
RBC antibody screen (4-cell)	4_Cell	Capture-R RS 4 Cells	Capture-R RS 4	
RBC antibody panel	Ab_ID	Capture-R RS Identification	Capture-R RS ID	
RBC antibody panel	ExtendDP	Capture-R RS Identification Extend1	Capture-R RS ID Extend1	

Assay Description	Assay Short Name	Generic name	Microplates used
	Capture-R Ready Screen		
RBC antibody panel	ExtendDN	Capture-R RS Identification Extend2	Capture-R RS ID Extend2
	Capture-CMV		
CMV detection	CMV	Capture-CMV	Capture-CMV

### Limitations, Warnings and Notes

#### lcon



**Limitation**: The purpose of the Monoclonal Control for the Weak\_D assays is to serve as a sample control. It is expected to indicate those sample related conditions that could lead to false positive well results and, therefore, a false positive interpretation for the test. When the Monoclonal Control well yields a positive result for a sample, then the instrument will give an interpretation of invalid.

Description



**Limitation:** The purpose of the Monoclonal Control for the Kell and Rh phenotyping assays is to serve as a sample control for Immucor low protein blood grouping reagents. It is expected to indicate those sample related conditions that could lead to spontaneous agglutination with low protein reagents and, therefore, a false positive interpretation for the test. When the Monoclonal Control well yields a positive result, then the instrument will give an interpretation of invalid.

#### lcon

### Description



**Limitation:** Certain blood sample factors have a high likelihood of causing no-typedetermined (NTD) ABO or Rh (D) interpretations. These include serological factors due to the inheritance of weakly expressed gene products by certain diseases (such as leukemia), through the transfusion or transplantation of ABO and Rh (D) allogeneic red blood cell products, or due to the patient's age. Instrument NTD results can also occur due to the interference caused by certain sample conditions, such as higher-than-normal levels of lipid, bilirubin, free plasma hemoglobin, clots, or aggregates.

It is possible that when two or more of these factors occur simultaneously, believable but erroneous test results (i.e., mistype) might occur. Forward only ABO-Rh testing has a higher risk of mistype due to the absence of the reverse type results. Hazardous mistypes may occur, such as an A sample being interpreted as group AB, or an Rh (D) negative sample being interpreted as Rh (D) positive. For this reason, ABO-Rh results should always be compared to the patient or donor's history.

Additionally, sample preparation and exclusion instructions contained in the relevant reagent direction circulars should be followed precisely. Before decisions about medical treatment are made or red blood cell products are released for transfusion, ABO and Rh (D) test results should be verified. This verification can consist of an immediate spin crossmatch or a result comparison to a second current or historical ABO and Rh (D) test by the same or an alternative method. This limitation applies to all ABO-Rh assays, both with and without a reverse test.



**Limitation:** The Reverse ABO Blood Group assay is intended to screen for ABO antibodies in plasma. This provides a presumptive indication of the red blood cell ABO group for the individual from which the plasma was collected. The absence of the corresponding forward group test subjects the results of this assay to several potential sources of error. These include serological factors due to the donor's age, immunological state, the presence of an IgM allo or autoantibody, transfusion or transplantation status or disease states. Sample related conditions such as higher-than-expected levels of lipids, bilirubin, free plasma hemoglobin, clots or aggregates may also be a factor. **This assay cannot be used (alone) to determine the ABO group for blood components intended for transfusion**.



**Limitation:** The red blood cell crossmatch assay (XMatch\_E/XMatch\_S) is intended only for the detection of incompatibilities due to IgG antibodies. The red blood cell crossmatch assay (XMatch\_E/XMatch\_S) is **not** intended for the detection of incompatibilities due to IgM antibodies, such as ABO incompatibilities.

#### Icon



Limitation: The NEO Iris software manages the results obtained based on Capture-R Ready-ID strip placement. To assure proper sample and well identification, it is necessary to load the Capture-R Ready-ID strips into the plate carrier, so that the left strip (L) (containing donors 1-8) is occupying an oddnumbered column position of the strip carrier. The strip must not be placed in the holder such that the left strip is occupying an even-numbered column position.

Description

#### **Correct Strip Positions**

**Column** Position

		1	2	3	4	5	6	7	8	9	10	11	12
	А	L	R	L	R	L	R	L	R	L	R	L	R
	В	L	R	L	R	L	R	L	R	L	R	L	R
	С	L	R	L	R	L	R	L	R	L	R	L	R
	D	L	R	L	R	L	R	L	R	L	R	L	R
	Е	L	R	L	R	L	R	L	R	L	R	L	R
ion	F	L	R	L	R	L	R	L	R	L	R	L	R
Posit	G	L	R	L	R	L	R	L	R	L	R	L	R
Row	Н	L	R	L	R	L	R	L	R	L	R	L	R

#### **Incorrect Strip Positions**

**Column Positions** 

		1	2	3	4	5	6	7	8	9	10	11	12
А	А		L	R	L	R	L	R	L	R	L	R	
	В		L	R	L	R	L	R	L	R	L	R	
	С		L	R	L	R	L	R	L	R	L	R	
	D		L	R	L	R	L	R	L	R	L	R	
E	E		L	R	L	R	L	R	L	R	L	R	
_	F		L	R	L	R	L	R	L	R	L	R	
ositior	G		L	R	L	R	L	R	L	R	L	R	
Row P	Н		L	R	L	R	L	R	L	R	L	R	

If this convention is not followed, the software will not recognize that the strips are misplaced and errors will occur. This scenario also applies to the Extend assays.

#### lcon

### Description



**Limitation:** Capture-R<sup>®</sup> Ready Indicator Red Cells should be used no more than 24 hours after a stir ball has been added to the vial. Vials of reagents, other than Indicator Red Cells, that have remained continuously on the NEO Iris for 72 hours (3 days) should be removed and replaced with fresh vials. Vials of reagents, other than Indicator Red Cells, that are removed from the NEO Iris when not in use and refrigerated can be used up to their expiration dates.



**Limitation:** In order to maintain the integrity of the DAT Positive Control Cell reagent after a stir ball has been added to the vial, Immucor recommends the following restrictions on reagent vial usage: (a) a maximum of 72 hours of continuous use on the NEO, or (b) a maximum of 7 days of intermittent use on the NEO Iris when the reagent vial is on-board daily for a maximum of 16 hours followed by intervening periods of refrigerated storage.



**Limitation:** The NEO Iris cannot reliably detect hemagglutination reactions that are graded as 1+ or less in test tube methodology.

**Warning:** When running the reagent quality control (QC) assays, only one lot of each reagent must be loaded onto the NEO.

**Warning:** The grading of reactions on the NEO Iris must only be regarded as an approximation when compared to off-line visual grading by laboratory technical staff.

**Warning:** For the QC3\_Cell assay, the correct position for the Capture-R<sup>®</sup> Ready-Screen<sup>®</sup> (3) strip in a Capture-R<sup>®</sup> Ready-Screen<sup>®</sup> (3) white plate frame is at strip position one (1). The instrument will not check for the correct placement of the Capture-R<sup>®</sup> Ready-Screen<sup>®</sup> (3) strip in the white plate frame. You must make sure of the correct placement of this strip yourself. Incorrect placement of this single strip can result in test material being pipetted onto the pipetting station surface, an overflow error in the washer, and a failed quality control. Unused Capture-R<sup>®</sup> Ready-Screen<sup>®</sup> (3) strips should not be present in any of the other positions of the plate. It is not required that strip positions two (2) to twelve (12) of the white plate frame be occupied by any strips at all. Those positions can be left empty, but all strip positions should be activated in the software. The *Strip Selection* tab in the *Plate Loading Tower* dialog must indicate that all twelve (12) strips are available even though only the strip at position (1) is present.



**Note:** The Ag\_K assay has a negative control well associated with each sample tested on the microplate. This control well contains the sample red blood cell suspension, Specimen Diluent and the Monoclonal Control reagent.

### lcon



**Note:** A negative control well containing Monoclonal Control is included with each sample tested using various assays to detect spontaneously agglutinated red blood cells and potential issues with the centrifugation and shake steps.

Description



**Note:** The full and forward-only blood grouping assays have been configured to require a strong positive (3-4+) reaction of the Anti-A, Anti-B, Anti-A,B and Anti-D reagents in order to return a positive well result. Weak positive reactivity with one of these reagents will return an equivocal "?" well result. This is so that clinical conditions resulting in 2+ reactivity or less can be made readily evident and may be investigated as needed according to your local on-site procedures and regulations. Such clinical conditions may include, but are not be limited to, ABO subgroups, D variants, and mixed-field reactions due to post bone marrow transplantation relapse or post blood transfusion.

# Assay Cutoffs

# List of Horizontal Assay Cutoffs

Horizontal Assays								
Assay Abbreviation	Reaction Grad		Lower Limit >	Upper Limit <=				
ABDCHECK		0	0	30				
		?	30	58				
	NOVACIONE Anti A	1+ (not reported)	N/A	N/A				
	NOVACIONE ANII-A	2+ (not reported)	N/A	N/A				
		3+ (not reported)	N/A	N/A				
		4+	58	100				
		0	0	30				
	NOVACLONE Anti-B	?	30	76				
	NOVACLONE Anti-D	1+(not reported)	N/A	N/A				
	NOVACLONE Diluent	2+(not reported)	N/A	N/A				
	Control	3+(not reported)	N/A	N/A				
		4+	76	100				
ABDCHECK2		0	0	30				
		?	30	58				
	ineme Clause Anti A	1+ (not reported)	N/A	N/A				
	Immucione Anti-A	2+ (not reported)	N/A	N/A				
		3+ (not reported)	N/A	N/A				
		4+	58	100				
		0	0	30				
	immuClone Anti-B	?	30	76				
	immuClone Rh-Hr Control	1+(not reported)	N/A	N/A				
		2+(not reported)	N/A	N/A				

Horizontal Assays						
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=		
		3+(not reported)	N/A	N/A		
		4+	76	100		
ABORH		0	0	30		
		?	30	58		
		1+ (not reported)	N/A	N/A		
	NOVACIONE ANII-A	2+ (not reported)	N/A	N/A		
		3+ (not reported)	N/A	N/A		
		4+	58	100		
		0	0	30		
	NOVACLONE Anti-B	?	30	76		
	NOVACLONE Anti-D	1+(not reported)	N/A	N/A		
	NOVACLONE Diluent	2+(not reported)	N/A	N/A		
	Control	3+(not reported)	N/A	N/A		
		4+	76	100		
		0	0	23		
		?	23	28		
	Referencells A1	1+	28	35		
	Referencells B	2+	35	50		
		3+	50	76		
		4+	76	100		
ABORH2		0	0	30		
		?	30	58		
		1+ (not reported)	N/A	N/A		
	Immucione Anti-A	2+ (not reported)	N/A	N/A		
		3+ (not reported)	N/A	N/A		
		4+	58	100		

	Horizontal Assays							
Assay	Reaction	Grade	Lower Limit >	Upper Limit <=				
Abbreviation								
		0	0	30				
	immuClone Anti-B	?	30	76				
	immuClone Anti-D rapid	1+(not reported)	N/A	N/A				
	immuClone Rh-Hr Control	2+(not reported)	N/A	N/A				
		3+(not reported)	N/A	N/A				
		4+	76	100				
		0	0	23				
		?	23	28				
	Referencells A1	1+	28	35				
	Referencells B	2+	35	50				
		3+	50	76				
		4+	76	100				
ABDLONG		0	0	30				
		?	30	58				
		1+ (not reported)	I Assays         Grade       Lower Limit >         0       30         ot reported)       N/A         23       28         24       35         50       50         76       0         28       35         50       76         0       30         not reported)       N/A         not reported)       N/A <td< td=""><td>N/A</td></td<>	N/A				
	NOVACLONE ANTI-A	2+ (not reported)	N/A	N/A				
		3+ (not reported)	N/A	N/A				
		4+	58	100				
	NOVACLONE Anti-B	0	0	30				
	NOVACLONE Anti-A,B	?	30	76				
	NOVACLONE Anti-D	1+(not reported)	N/A	N/A				
	immuClone Anti-D rapid	2+(not reported)	N/A	N/A				
	immuClone Rh-Hr Control	3+(not reported)	N/A	N/A				
		4+	76	100				
	Referencells A1	0	0	23				
	Referencells A2	?	23	28				
Horizontal Assays								
-----------------------	---	-------------------	---------------	----------------				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=				
	Referencells B	1+	28	35				
	Referencells O Autocontrol	2+	35	50				
		3+	50	76				
		4+	76	100				
ABDLONG2		0	0	30				
		?	30	58				
		1+ (not reported)	N/A	N/A				
	ImmuClone Anti-A	2+ (not reported)	N/A	N/A				
		3+ (not reported)	N/A	N/A				
		4+	58	100				
	immuClone Anti-B immuClone Anti-A,B NOVACLONE Anti-D immuClone Anti-D rapid immuClone Rh-Hr Control	0	0	30				
		?	30	76				
		1+(not reported)	N/A	N/A				
		2+(not reported)	N/A	N/A				
		3+(not reported)	N/A	N/A				
		4+	76	100				
	Referencells A1 Referencells A2 Referencells B	0	0	23				
		?	23	28				
		1+	28	35				
	Referencells O	2+	35	50				
	Autocontrol	3+	50	76				
		4+	76	100				
BABY_BG		0	0	30				
	immuClone Anti-A	?	30	58				
	NOVACLONE Anti-A	1+ (not reported)	N/A	N/A				
		2+ (not reported)	N/A	N/A				

Horizontal Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		3+ (not reported)	N/A	N/A
		4+	58	100
	immuClone Anti-B	0	0	30
	NOVACLONE Anti-B	?	30	76
	immuClone Anti-A,B	1+(not reported)	N/A	N/A
	NOVACLONE Anti-A,B	2+(not reported)	N/A	N/A
	NOVACLONE Anti-D	3+(not reported)	N/A	N/A
	immuClone Anti-D rapid immuClone Rh-Hr Control	4+	76	100
A_subg	A_subg Anti-A1 (Lectin) immuClone Anti-H immuClone Rh-Hr Control	0	0	23
		?	23	28
		1+	28	35
		2+	35	50
		3+	50	80
		4+	80	100
AB_CTR		0	0	30
		?	30	58
	NOVACIONE Anti A	1+ (not reported)	N/A	N/A
	NOVACLONE ANTI-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
	NOVACLONE Anti-B	?	30	76
	Control	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A

Horizontal Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		4+	76	100
AB_CTR2		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	Immucione Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
		?	30	76
	immuClone Anti-B immuClone Rh-Hr Control	1+ (not reported)	N/A	N/A
ır		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
RHFORMEL		0	0	23
	immuClone (1) Anti-C	?	23	75
	immuClone (2) Anti-C	1+ (not reported)	N/A	N/A
	immuClone (1) Anti-c immuClone (2) Anti-c	2+ (not reported)	N/A	N/A
		3+	75	80
immuClone (2) Anti-E immuClone (1) Anti-e immuClone (2) Anti-e immuClone (2) Anti-e	4+	80	100	
		0	0	23
	Immucione (1) Anti-K	?	23	35
		1+ (not reported)	N/A	N/A

Horizontal Assays				
Assay	Reaction	Grade	Lower Limit >	Upper Limit <=
Abbreviation				
		2+	35	50
		3+	50	80
		4+	80	100
		0	0	23
		?	23	50
	Automated immuClone	1+ (not reported)	N/A	N/A
	Anti-K	2+ (not reported)	N/A	N/A
		3+	50	80
		4+	80	100
PHENO16	immuClone (1) Anti-C	0	0	23
	immuClone (1) Anti-c	?	23	75
	immuClone (1) Anti-E	1+ (not reported)	N/A	N/A
	immuClone (1) Anti-e	2+ (not reported)	N/A	N/A
	immuClone Rh-Hr Control	3+	75	80
		4+	80	100
	immuClone (1) Anti-K	0	0	23
		?	23	35
		1+ (not reported)	N/A	N/A
		2+	35	50
		3+	50	80
		4+	80	100
AG_K		0	0	23
		?	23	35
	ImmuClone (1) Anti-K	1+ (not reported)	N/A	N/A
		2+	35	50
		3+	50	80

Horizontal Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		4+	80	100
		0	0	23
		?	23	75
	immuClone Rh-Hr Control	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	100
Kell		0	0	23
	immuClone (1) Anti-K	?	23	35
		1+ (not reported)	N/A	N/A
		2+	35	50
		3+	50	80
		4+	80	100
	Automated immuClone Anti-K	0	0	23
		?	23	50
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	50	80
		4+	80	100
		0	0	23
		?	23	75
	immuClone Rh-Hr Control	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	100
AG_M	immuClone Anti-M	0	0	23

Horizontal Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		?	23	75
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	75	100
	immuClone Rh-Hr Control	0	0	23
		?	23	80
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	80	100
Cw_lgM	immuClone Anti-Cw	0	0	23
		?	23	75
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	100
		0	0	23
		?	23	80
	immuClone Rh-Hr Control	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	80	100

## List of Vertical Assay Cutoffs

Vertical Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
ABDCHECK		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	NOVACLONE ANTI-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
	NOVACLONE Anti-B	?	30	76
	NOVACLONE Anti-D NOVACLONE Diluent Control	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+(not reported)	N/A	N/A
		4+	76	100
ABDCHECK2		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	ImmuCione Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
	immuClone Anti-B	?	30	76
	immuClone Anti-D rapid	1+(not reported)	N/A	N/A
	immuClone Rh-Hr	2+(not reported)	N/A	N/A
	Control	3+(not reported)	N/A	N/A
		4+	76	100

Vertical Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
ABDFULL		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	NOVACLONE ANTI-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
	NOVACIONE Anti-B	0	0	30
	NOVACLONE Anti-A,B	?	30	76
	NOVACLONE Anti-D immuClone Anti-D rapid NOVACLONE Diluent Control	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+(not reported)	N/A	N/A
		4+	76	100
	Referencells A1	0	0	23
		?	23	28
	Referencells B	1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
ABDFULL2		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
	immuClone Anti-B	0	0	30
	immuClone Anti-A,B	?	30	76

Vertical Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
	immuClone Anti-D rapid	1+(not reported)	N/A	N/A
	immuClone Rh-Hr	2+(not reported)	N/A	N/A
	Control	3+(not reported)	N/A	N/A
		4+	76	100
		0	0	23
		?	23	28
	Referencells A1	1+	28	35
	Referencells B	2+	35	50
		3+	50	76
		4+	76	100
ABOD12_I	immuClone Anti-A	0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
	immuClone Anti-B	?	30	76
	immuClone Anti-A,B	1+(not reported)	N/A	N/A
	immuClone Anti-D rapid	2+(not reported)	N/A	N/A
	Control	3+(not reported)	N/A	N/A
		4+	76	100
		0	0	30
	Empty Row	? (not reported)	N/A	N/A
		Positive	30	100
	Referencells A1	0	0	23

Vertical Assays				
Assay	Reaction	Grade	Lower Limit >	Upper Limit <=
Abbreviation				
	(GRA1)	?	23	28
	Referencells B	1+	28	35
	(GRB)	2+	35	50
		3+	50	76
		4+	76	100
ABOD12NC		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	NOVACLONE ANTI-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
	NOVACLONE Anti-B	0	0	30
	NOVACLONE Anti-A,B NOVACLONE Anti-D NOVACLONE Diluent Control	?	30	76
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+(not reported)	N/A	N/A
		4+	76	100
		0	0	30
		? (not reported)	N/A	N/A
		Positive	30	100
		0	0	23
	Referencells A1	?	23	28
	(GRA1)	1+	28	35
	Referencells B	2+	35	50
	(GRB)	3+	50	76
		4+	76	100

Vertical Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
Baby_BG		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	NOVACLONE ANTI-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
	NOVACIONE Anti P	0	0	30
	NOVACLONE Anti-A,B	?	30	76
	NOVACLONE Anti-D immuClone Anti-D rapid NOVACLONE Diluent Control	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+(not reported)	N/A	N/A
		4+	76	100
Baby_BG2		0	0	30
		?	30	58
	immuClana Anti A	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
	immuClone Anti-B	?	30	76
	immuClone Anti-A,B	1+(not reported)	N/A	N/A
	immuClone Anti-D rapid	2+(not reported)	N/A	N/A
	immuClone Rh-Hr	3+(not reported)	N/A	N/A
		4+	76	100
ABODBB_I		0	0	30
	immuClone Anti-A	?	30	58

Vertical Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
	immuClone Anti-B	?	30	76
	immuClone Anti-A,B	1+(not reported)	N/A	N/A
	immuClone Anti-D rapid	2+(not reported)	N/A	N/A
	Control	3+(not reported)	N/A	N/A
		4+	76	100
A_subg		0	0	23
	Anti-A1 (Lectin) immuClone Anti-H immuClone Rh-Hr Control	?	23	28
		1+	28	35
		2+	35	50
		3+	50	80
		4+	80	100
AG_H		0	0	23
		?	23	28
	immuClone Anti-H	1+	28	35
	ImmuCione Rh-Hr	2+	35	50
		3+	50	80
		4+	80	100
Pheno12	immuClone (1) Anti-C	0	0	23
	immuClone (1) Anti-c	?	23	75
	immuClone (1) Anti-E	1+ (not reported)	N/A	N/A
	immuClone (1) Anti-e	2+ (not reported)	N/A	N/A

Vertical Assays					
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=	
	immuClone Rh-Hr	3+	75	80	
	Control	4+	80	100	
		0	0	23	
		?	23	35	
	immuClone (1) Anti-K	1+ (not reported)	N/A	N/A	
		2+	35	50	
		3+	50	80	
		4+	80	100	
Pheno	immuClone (1) Anti-C	0	0	23	
	immuClone (1) Anti-c immuClone (1) Anti-E immuClone (1) Anti-e immuClone Rh-Hr Control	?	23	75	
		1+ (not reported)	N/A	N/A	
		2+ (not reported)	N/A	N/A	
		3+	75	80	
		4+	80	100	
		0	0	23	
		?	23	50	
	Automated immuClone	1+ (not reported)	N/A	N/A	
	Anti-K	2+ (not reported)	N/A	N/A	
		3+	50	80	
		4+	80	100	
AG_K		0	0	23	
		?	23	50	
	Automated immuClone	1+(not reported)	N/A	N/A	
		2+(not reported)	N/A	N/A	
		3+	50	80	

Vertical Assays					
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=	
		4+	80	100	
		0	0	23	
		?	23	50	
	immuClone Rh-Hr	1+ (not reported)	N/A	N/A	
	Control	2+ (not reported)	N/A	N/A	
		3+	50	80	
		4+	80	100	
AG_M		0	0	23	
	immuClone Anti-M	?	23	75	
		1+ (not reported)	N/A	N/A	
		2+ (not reported)	N/A	N/A	
		3+ (not reported)	N/A	N/A	
		4+	75	100	
		0	0	23	
		?	23	80	
	immuClone Rh-Hr	1+ (not reported)	N/A	N/A	
	Control	2+ (not reported)	N/A	N/A	
		3+ (not reported)	N/A	N/A	
		4+	80	100	
Cw_lgM		0	0	23	
		?	23	75	
		1+ (not reported)	N/A	N/A	
	Immucione Anti-Cw	2+ (not reported)	N/A	N/A	
		3+	75	80	
		4+	80	100	
	immuClone Rh-Hr	0	0	23	

Vertical Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
	Control	?	23	80
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	80	100
Pool_Cell		0	0	20
2_Cell		?	20	30
4_Cell		1+	30	45
	All	2+	45	65
		3+	65	90
		4+	90	100
3_Cell		0	0	20
		?	20	30
	Screening wells	1+	30	45
		2+	45	65
		3+	65	90
		4+	90	100
		0	0	30
		?	-	-
	IgG coated positive	1+	30	45
	control well	2+	45	65
		3+	65	90
		4+	90	100
Ab_ID		0	0	20
ExtendDP	Sample	?	20	30
ExtendDN		1+	30	45

Vertical Assays					
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=	
		2+	45	65	
		3+	65	90	
		4+	90	100	
		0	0	30	
		?	-	-	
		1+	30	45	
	Positive Control	2+	45	65	
		3+	65	90	
		4+	90	100	
		0	0	30	
		?	-	-	
	Negative Control	1+	30	45	
		2+	45	65	
		3+	65	90	
		4+	90	100	
DAT		0	0	20	
		?	20	40	
		1+	40	50	
	All	2+	50	72	
		3+	72	90	
		4+	90	100	
XMATCH_E		0	0	5	
XMATCH_S		?	5	15	
	All	1+	15	30	
		2+	30	60	
		3+	60	90	

Vertical Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		4+	90	100
Weak_D1,		0	0	20
Weak_D2		?	20	40
		1+	40	50
	Anti-D	2+	50	72
		3+	72	90
		4+	90	100
		0	0	20
		?	20	40
	immuClone Rh-Hr Control	1+	40	50
		2+	50	72
		3+	72	90
		4+	90	100
AG_Cw	Anti-Cw micro	0	0	20
		?	20	40
		1+	40	50
		2+	50	72
		3+	72	90
		4+	90	100
		0	0	20
		?	20	40
		1+	40	50
	Negative Control micro	2+	50	72
		3+	72	90
		4+	90	100
AG_Fya	Anti-Fy(a) micro	0	0	20

Vertical Assays					
Assay	Reaction	Grade	Lower Limit >	Upper Limit <=	
Abbreviation					
		?	20	40	
		1+	40	50	
		2+	50	72	
		3+	72	90	
		4+	90	100	
		0	0	20	
		?	20	40	
		1+	40	50	
	Negative Control micro	2+	50	72	
		3+	72	90	
		4+	90	100	
AG_Fyb		0	0	20	
		?	20	40	
		1+	40	50	
	Anti-Fy(b) micro	2+	50	72	
		3+	72	90	
		4+	90	100	
		0	0	20	
		?	20	40	
		1+	40	50	
	Negative Control micro	2+	50	72	
		3+	72	90	
		4+	90	100	
AG_Jka		0	0	20	
	Anti-Jk(a) micro	?	20	40	
		1+	40	50	

Vertical Assays				
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		2+	50	72
		3+	72	90
		4+	90	100
		0	0	20
		?	20	40
	No softing Company and and	1+	40	50
	Negative Control micro	2+	50	72
		3+	72	90
		4+	90	100
AG_Jkb		0	0	20
		?	20	40
	Anti-Jk(b) micro	1+	40	50
		2+	50	72
		3+	72	90
		4+	90	100
		0	0	20
		?	20	40
		1+	40	50
	Negative Control micro	2+	50	72
		3+	72	90
		4+	90	100
AG_k_li		0	0	20
		?	20	40
	Anti k (Collono) micro	1+	40	50
		2+	50	72
		3+	72	90
		4+	90	100

Vertical Assays					
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=	
		0	0	20	
		?	20	40	
		1+	40	50	
	Negative Control micro	2+	50	72	
		3+	72	90	
		4+	90	100	
AG_S		0	0	20	
		?	20	40	
		1+	40	50	
	Anti-S micro	2+	50	72	
		3+	72	90	
		4+	90	100	
	Negative Control micro	0	0	20	
		?	20	40	
		1+	40	50	
		2+	50	72	
		3+	72	90	
		4+	90	100	
AG_s_li		0	0	20	
		?	20	40	
	Anti a miana	1+	40	50	
	Anti-s micro	2+	50	72	
		3+	72	90	
		4+	90	100	
		0	0	20	
	Negative Control micro	?	20	40	
		1+	40	50	

Vertical Assays					
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=	
		2+	50	72	
		3+	72	90	
		4+	90	100	
ANTIGENS	Anti-Fy(a) micro	0	0	20	
	Anti-Fy(b) micro	?	20	40	
	Anti-Jk(a) micro	1+	40	50	
	Anti-Jk(b) micro	2+	50	72	
	Anti-S micro Anti-s micro	3+	72	90	
		4+	90	100	
		0	0	20	
		?	20	40	
Negative Control mic		1+	40	50	
	Negative Control micro	2+	50	72	
		3+	72	90	
		4+	90	100	

N/A = Not Applicable

Assay Abbreviation	Reaction	Neg <=	Pos >	
CMV	Not Applicable	45.5	45.5	

## Assay Reagent Component Grid

Reagents & Microplates	Pool_Cell	2_Cell	3_Cell	4_Cell	Ab_ID	ExtendDP	ExtendDN	CMV
Capture LISS	Х	Х	х	х	х	х	х	х
Capture-R Ready Indicator Red Cells	х	Х	х	Х	х	х	х	
Capture-R Positive Control Serum (Weak)	х	х		х	х	х	х	
Capture-R Negative Control Serum	х	Х		Х	Х	Х	Х	
Capture-R Ready-Screen (Pooled Cells)	х							
Capture-R Ready-Screen(I and II)		Х						
Capture-R Ready-Screen (3)			х					
Capture-R Ready-Screen (4)				Х				
Capture-R Ready-ID					Х			
Capture-R Ready-ID Extend I						Х		
Capture-R Ready-ID Extend II							Х	
Capture-CMV Indicator Red Cells								х
Capture-CMV Positive Control Serum (Weak)								x
Capture-CMV Negative Control Serum								х
Capture-CMV								х

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Capture-R Select		
Direct Antiglobulin Test (DAT)	DAT	N/A (only controls)	Capture-R Select Plates
Weak D testing	Weak_D1	1. NOVACLONE Anti-D	Capture-R Select Plates
Weak D testing	Weak_D2	<ol> <li>NOVACLONE Anti-D</li> <li>immuClone Rh-Hr Control</li> </ol>	Capture-R Select Plates
RBC crossmatch for donor-samples	XMatch_E	N/A (only controls)	Capture-R Select Plates
in EDTA tubes			
RBC crossmatch for segment donor-samples	XMatch_S	N/A (only controls)	Capture-R Select Plates

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Capture Select - Rare Antigens		
Capture-R Cw Antigen	AG_Cw	<ol> <li>Anti-Cw micro</li> <li>Negative Control micro</li> </ol>	Capture-R Select Plates
Capture-R Fya Antigen	AG_Fya	<ol> <li>Anti-Fya micro</li> <li>Negative Control micro</li> </ol>	Capture-R Select Plates
Capture-R Fyb Antigen	AG_Fyb	<ol> <li>Anti-Fyb micro</li> <li>Negative Control micro</li> </ol>	Capture-R Select Plates
Capture-R Jka Antigen	AG_Jka	<ol> <li>Anti-Jka</li> <li>Negative Control micro</li> </ol>	Capture-R Select Plates
Capture-R Jkb Antigen	AG_Jkb	<ol> <li>Anti-Jkb micro</li> <li>Negative Control micro</li> </ol>	Capture-R Select Plates

Assay Description	Assay Short Name	Used Reagents Microplates used
	Capture Select - Rare Antigens	
Capture-R k Antigen	AG_k_li	1. Anti-Cellano microCapture-R Select Plates2. Negative Control micro
Capture-R S Antigen	AG_S	1. Anti-S micro     Capture-R Select Plates       2. Negative Control micro
Capture-R s Antigen	AG_s_li	1. Anti-s micro     Capture-R Select Plates       2. Negative Control micro
Capture-R Antigens	ANTIGENS	<ol> <li>Anti-Fya micro</li> <li>Anti-Fyb micro</li> <li>Anti-Jka micro</li> <li>Anti-Jkb micro</li> <li>Anti-Jkb micro</li> <li>Anti-S micro</li> <li>Anti-s micro</li> <li>Empty</li> <li>Negative Control micro</li> </ol>

## Assay Procedural Steps

## **Before You Begin**



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button Abbreviation		Brief Synopsis of Assay Procedural Steps
ABORH	1.	Bring all reagents and blood samples to 18–30°C before testing.
ABORH_2	2.	Centrifuge the blood samples to separate the plasma from the red blood cells. Remove the caps from the blood sample tubes.
	3.	Select the necessary number of untreated microplates for use.
	4.	Remove reagent vial caps.
	5.	Add one stirball to each new vial of Referencells <sup>®</sup> A1 and B to be used. Gently agitate each vial to resuspend the red blood cells.
	6.	Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7.	Assign the ABORH or the ABORH_2 assay to the blood samples, either manually or following the upload worklist procedure.
	8.	Start the ABORH or the ABORH_2 assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the ABORH or the ABORH_2 assay, and records and interprets blood sample results.
	9.	At the completion of the NEO Iris ABORH or the ABORH_2 assay, click the Results button on the main menu bar to access the blood sample results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
ABDLONG	1. Bring all reagents and blood samples to 18–30°C before testing.
ABDLONG2	2. Centrifuge the blood samples to separate the plasma from the red blood cells. Remove the caps from the blood sample tubes.
	3. Select the necessary number of untreated microplates for use.
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Referencells <sup>®</sup> A1, A2, B and O to be used. Gently agitate each vial to resuspend the red blood cells.
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the ABDLONG or the ABDLONG2 assay to the blood samples, either manually or following the upload worklist procedure.
	<ol> <li>Start the ABDLONG or the ABDLONG2 assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the ABDLONG or the ABDLONG2 assay, and records and interprets blood sample results.</li> </ol>
	9. At the completion of the NEO Iris ABDLONG or the ABDLONG2 assay, click the Results button on the main menu bar to access the blood sample results

Assay Button Abbreviation		Brief Synopsis of Assay Procedural Steps
ABDCHECK ABDCHECK2	1.	Bring all reagents, donor segment blood samples and blood samples (if applicable) to 18–30°C before testing.
	2.	Centrifuge the blood sample collection tubes to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge these samples.
	3.	Select the necessary number of untreated microplates for use.
	4.	Remove reagent vial caps.
	5.	Load reagents, microplates, and blood samples (and/or donor samples) onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	6.	Assign the ABDCHECK or the ABDCHECK2 assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	7.	Start the ABDCHECK or the ABDCHECK2 assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the ABDCHECK or the ABDCHECK2 assay, and records and interprets the blood sample (and/or donor samples) results.
	8.	At the completion of the NEO Iris ABDCHECK or the ABDCHECK2 assay, click the Results button on the main menu bar to access the blood sample (and/or donor samples) results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
BABY_BG	<ol> <li>Bring all reagents, donor segment blood samples and blood samples (if applicable) to 18–30°C before testing.</li> </ol>
	2. Centrifuge the blood sample collection tubes to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge these samples.
	3. Select the necessary number of untreated microplates for use.
	4. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples) onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	6. Assign the BABY_BG assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	7. Start the BABY_BG assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the BABY_BG assay, and records and interprets the blood sample (and/or donor samples) results.
	8. At the completion of the NEO Iris BABY_BG assay, click the Results button on the main menu bar to access the blood sample (and/or donor samples) results.

Assay Button Abbreviation		Brief Synopsis of Assay Procedural Steps
AB_CTR AB_CTR2	1.	Bring all reagents, donor segment blood samples and blood samples (if applicable) to 18–30°C before testing.
	2.	Centrifuge the blood sample collection tubes to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge these samples.
	3.	Select the necessary number of untreated microplates for use.
	4.	Remove reagent vial caps.
	5.	Load reagents, microplates, and blood samples (and/or donor samples) onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	6.	Assign the AB_CTR or the AB_CTR2 assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	7.	Start the AB_CTR or the AB_CTR2 assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the AB_CTR or the AB_CTR2 assay, and records and interprets the blood sample (and/or donor samples) results.
	8.	At the completion of the NEO Iris AB_CTR or the AB_CTR2 assay, click the Results button on the main menu bar to access the blood sample (and/or donor samples) results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
A_subg	<ol> <li>Bring all reagents, donor segment blood samples and blood samples (if applicable) to 18–30°C before testing.</li> </ol>
	2. Centrifuge the blood sample collection tubes to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge these samples.
	3. Select the necessary number of untreated microplates for use.
	4. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples) onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	<ol> <li>Assign the A_subg assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.</li> </ol>
	<ol> <li>Start the A_subg assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the A_subg assay, and records and interprets the blood sample (and/or donor samples) results.</li> </ol>
	8. At the completion of the NEO Iris A_subg assay, click the Results button on the main menu bar to access the blood sample (and/or donor samples) results.

Assay Button Abbreviation		Brief Synopsis of Assay Procedural Steps
RHFORMEL	1.	Bring all reagents and blood samples to 18–30°C before testing.
PHENO16	2.	Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.
	3.	Remove reagent vial caps.
	4.	Load reagents, microplates, and blood samples (and/or donor samples), onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	5.	Assign the RHFORMEL or the PHENO16 assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	6.	Start the RHFORMEL or the PHENO16 assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the RHFORMEL or the PHENO16 assay, and records and interprets blood samples results.
	7.	At the completion of the NEO Iris RHFORMEL or the PHENO16 assay, click the Results button on the main menu bar to access the blood samples results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
AG_K	1. Bring all reagents and blood samples to 18–30°C before testing.
Kell	2. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.
	3. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples), onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	5. Assign the AG_K or the Kell assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	<ol> <li>Start the AG_K or the Kell assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the AG_K or the Kell assay, and records and interprets plasma/serum sample results.</li> </ol>
	7. At the completion of the NEO Iris AG_K or the Kell assay, click the Results button on the main menu bar to access the blood samples results.
AG_M	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.
	3. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples), onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	5. Assign the AG_M assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	<ol> <li>Start the AG_M assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the AG_M assay, and records and interprets blood samples results.</li> </ol>
	7. At the completion of the NEO Iris AG_M assay, click the Results button on the main menu bar to access the blood samples results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
Cw_lgM	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.
	3. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples), onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	5. Assign the Cw assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	6. Start the Cw assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the Cw assay, and records and interprets blood samples results.
	7. At the completion of the NEO Iris Cw assay, click the Results button on the main menu bar to access the blood samples results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
ABDCHECK ABDCHECK2	1. Bring all reagents, donor segment blood samples and blood samples (if applicable) to 18–30°C before testing.
	2. Centrifuge the blood sample collection tubes to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge these samples.
	3. Select the necessary number of untreated microplates for use.
	4. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples) onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	6. Assign the ABDCHECK or the ABDCHECK2 assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	<ol> <li>Start the ABDCHECK or the ABDCHECK2 assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the ABDCHECK or the ABDCHECK2 assay, and records and interprets the blood sample (and/or donor samples) results.</li> </ol>
	8. At the completion of the NEO Iris ABDCHECK or the ABDCHECK2 assay, click the Results button on the main menu bar to access the blood sample (and/or donor samples) results.

Assay Button Abbreviation		Brief Synopsis of Assay Procedural Steps
ABDFULL	1.	Bring all reagents and blood samples to 18–30°C before testing.
ABDFULL2 ABOD12_I	2.	Centrifuge the blood samples to separate the plasma from the red blood cells. Remove the caps from the blood sample tubes.
ABOD12_NC	3.	Select the blood samples to separate the plasma from the red blood cells. Remove the caps from the blood sample tubes.
	4.	Remove reagent vial caps.
	5.	Add one stirball to each new vial of Referencells® A1 and B to be used. Gently agitate each vial to resuspend the red blood cells.
	6.	Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7.	Assign the ABFULL, ABDFULL2, ABOD12_I or the ABOD12_NC assay to the blood samples, either manually or following the upload worklist procedure.
	8.	Start the ABFULL, ABDFULL2, ABOD12_I or the ABOD12_NC assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the ABFULL, ABDFULL2, ABOD12_I or the ABOD12_NC assay, and records and interprets blood sample results.
	9.	At the completion of the NEO Iris ABFULL, ABDFULL2, ABOD12_I or the ABOD12_NC assay, click the Results button on the main menu bar to access the blood sample results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps	
Baby_BG       1.         Baby_BG2       2.         ABODBB_I       3.         4.       5.         6.       7.         8.       8.	1. Bring all reagents, donor segment blood samples and blood samples (if applicable) to 18–30°C before testing.	
	2. Centrifuge the blood sample collection tubes to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge these samples.	
	3. Select the necessary number of untreated microplates for use.	
	4. Remove reagent vial caps.	
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples) onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>	
	<ol> <li>Assign the Baby_BG, Baby_BG2 or the ABODBB_I assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.</li> </ol>	
	7. Start the Baby_BG, Baby_BG2 or the ABODBB_I assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the Baby_BG, Baby_BG2 or the ABODBB_I assay, and records and interprets the blood sample (and/or donor samples) results.	
	8. At the completion of the NEO Iris Baby_BG, Baby_BG2 or the ABODBB_I assay, click the Results button on the main menu bar to access the blood sample (and/or donor samples) results.	
Assay Button Abbreviation		Brief Synopsis of Assay Procedural Steps
------------------------------	----	---
PHENO12	1.	Bring all reagents and blood samples to 18–30°C before testing.
PHENO	2.	Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.
	3.	Remove reagent vial caps.
	4.	Load reagents, microplates, and blood samples (and/or donor samples), onto the NEO Iris.
	5.	Follow the procedures in Chapter 6 – Instrument Testing Operation.
	6.	Assign the PHENO12 or the PHENO assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	7.	Start the PHENO12 or the PHENO assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the PHENO12 or the PHENO assay, and records and interprets blood samples results.
	8.	At the completion of the NEO Iris PHENO12 or the PHENO assay, click the Results button on the main menu bar to access the blood samples results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
AG_K	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.
	3. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples), onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	5. Assign the AG_K assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	<ol> <li>Start the AG_K assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the AG_K assay, and records and interprets blood samples results.</li> </ol>
	7. At the completion of the NEO Iris AG_K assay, click the Results button on the main menu bar to access the blood samples results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
A_subg	<ol> <li>Bring all reagents, donor segment blood samples and blood samples (if applicable) to 18–30°C before testing.</li> </ol>
	2. Centrifuge the blood sample collection tubes to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge these samples.
	3. Select the necessary number of untreated microplates for use.
	4. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples) onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	6. Assign the A_subg assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	<ol> <li>Start the A_subg assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the A_subg assay, and records and interprets the blood sample (and/or donor samples) results.</li> </ol>
	8. At the completion of the NEO Iris A_subg assay, click the Results button on the main menu bar to access the blood sample (and/or donor samples) results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
AG_H	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.
	3. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples), onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	5. Assign the AG_H assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	<ol> <li>Start the AG_H assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the AG_H assay, and records and interprets blood samples results.</li> </ol>
	7. At the completion of the NEO Iris AG_H assay, click the Results button on the main menu bar to access the blood samples results.
AG_M	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.
	3. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples), onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	<ol> <li>Assign the AG_M assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.</li> </ol>
	<ol> <li>Start the AG_M assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the AG_M assay, and records and interprets blood samples results.</li> </ol>
	7. At the completion of the NEO Iris AG_M assay, click the Results button on the main menu bar to access the blood samples results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
Cw_lgM	. Bring all reagents and blood samples to 18–30°C before testing.
	. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.
	. Remove reagent vial caps.
	. Load reagents, microplates, and blood samples (and/or donor samples), onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	. Assign the Cw assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	. Start the Cw assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the Cw assay, and records and interprets blood samples results.
	. At the completion of the NEO Iris Cw assay, click the Results button on the main menu bar to access the blood samples results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
Pool_Cell	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the blood samples from the red blood cells/clot. Remove the caps from the blood sample tubes.
	3. Remove the Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (Pooled Cells) microplate frame and the desired number of Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (Pooled Cells) strips from the pouch.
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the Pool_Cell assay to the blood samples, either manually or following the upload worklist procedure.
	<ol> <li>Start the Pool_Cell assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the Pool_Cell assay, and records and interprets the blood sample results.</li> </ol>
	9. At the completion of the NEO Iris Pool_Cell assay, click the Results button on the main menu bar to access the blood sample results.

Assay Button Abbreviation		Brief Synopsis of Assay Procedural Steps
2_Cell	1.	Bring all reagents and blood samples to 18–30°C before testing.
	2.	Centrifuge the blood samples to separate the blood samples from the red blood cells/clot. Remove the caps from the blood sample tubes.
	3.	Remove the Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (I and II) microplate frame and the desired number of Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (I and II) strips from the pouch.
	4.	Remove reagent vial caps.
	5.	Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	6.	Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7.	Assign the 2_Cell assay to the blood samples, either manually or following the upload worklist procedure.
	8.	Start the 2_Cell assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the 2_Cell assay, and records and interprets the blood sample results.
	9.	At the completion of the NEO Iris 2_Cell assay, click the Results button on the main menu bar to access the blood sample results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
3_Cell	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the blood samples from the red blood cells/clot. Remove the caps from the blood sample tubes.
	3. Remove the Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (3) microplate frame and the desired number of Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (3) strips from the pouch.
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the 3_Cell assay to the blood samples, either manually or following the upload worklist procedure.
	<ol> <li>Start the 3_Cell assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the 3_Cell assay, and records and interprets the blood sample results.</li> </ol>
	9. At the completion of the NEO Iris 3_Cell assay, click the Results button on the main menu bar to access the blood sample results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
QC3_Cell	1. Bring all reagents to 18–30°C before testing.
	2. Remove the Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (3) plate frame and one (1) Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (3) strip from the pouch. Place the strip in position one (1).
	3. Remove reagent vial caps.
	4. Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	<ol> <li>Load reagents and plates onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	<ol> <li>Start the QC3_Cell assay following the procedures in Chapter 10 – Maintaining the NEO. The NEO Iris automatically performs the QC3_Cell assay, and records and interprets the QC3_Cell results.</li> </ol>
	<ol> <li>At the completion of the NEO Iris QC3_Cell assay, press the Results button on the main menu bar to access the results. A QC3_Cell report will print automatically.</li> </ol>

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
4_Cell	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the blood samples from the red blood cells/clot. Remove the caps from the blood sample tubes.
	3. Remove the Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (4) microplate frame and the desired number of Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (4) strips from the pouch.
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the 4_Cell assay to the blood samples, either manually or following the upload worklist procedure.
	<ol> <li>Start the 4_Cell assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the 4_Cell assay, and records and interprets the blood sample results.</li> </ol>
	9. At the completion of the NEO Iris 4_Cell assay, click the Results button on the main menu bar to access the blood sample results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
Ab_ID	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the blood samples from the red blood cells/clot. Remove the caps from the blood sample tubes.
	3. Remove the Capture-R <sup>®</sup> Ready-ID <sup>®</sup> microplate frame and the desired number of Capture-R <sup>®</sup> Ready-ID <sup>®</sup> strips from the pouch.
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the Ab_ID assay to the blood samples, either manually or following the upload worklist procedure.
	<ol> <li>Start the Ab_ID assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the Ab_ID assay, and records the blood sample results.</li> </ol>
	9. At the completion of the NEO Iris Ab_ID assay, click the Results button on the main menu bar to access the blood sample results. If any wells are reported positive, you must cross-reference this data with the lot-specific Capture-R <sup>®</sup> Ready-ID <sup>®</sup> Master List to determine the antibody identification (if any exists).

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
ExtendDP	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the blood samples from the red blood cells/clot. Remove the caps from the blood sample tubes.
	3. Remove the Capture-R <sup>®</sup> Ready-ID <sup>®</sup> Extend I microplate frame and the desired number of Capture-R <sup>®</sup> Ready-ID <sup>®</sup> Extend I strips from the pouch.
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the ExtendDP assay to the blood samples, either manually or following the upload worklist procedure.
	8. Start the ExtendDP assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the ExtendDP assay, and records the blood sample results.
	9. At the completion of the NEO Iris ExtendDP assay, click the Results button on the main menu bar to access the blood sample results. If any wells are reported positive, you must cross-reference this data with the lot-specific Capture-R <sup>®</sup> Ready-ID <sup>®</sup> Extend I Master List to determine the antibody identification (if any exists).

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
ExtendDN	1. Bring all reagents and blood samples to 18–30°C before testing.
	<ol> <li>Centrifuge the blood samples to separate the blood samples from the red blood cells/clot. Remove the caps from the blood sample tubes.</li> </ol>
	<ol> <li>Remove the Capture-R<sup>®</sup> Ready-ID<sup>®</sup> Extend II microplate frame and the desired number of Capture-R<sup>®</sup> Ready-ID<sup>®</sup> Extend II strips from the pouch.</li> </ol>
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	5. Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the ExtendDN assay to the blood samples, either manually or following the upload worklist procedure.
	<ol> <li>Start the ExtendDN assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the ExtendDN assay, and records the blood sample results.</li> </ol>
	At the completion of the NEO Iris ExtendDN assay, click the Results button on the main menu bar to access the blood sample results. If any wells are reported positive, you must cross-reference this data with the lot-specific Capture-R <sup>®</sup> Ready-ID <sup>®</sup> Extend II Master List to determine the antibody identification (if any exists).

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
DAT	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood sample collection tubes to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge these samples.
	3. Remove the Capture-R <sup>®</sup> Select microplate frame and the desired number of Capture-R <sup>®</sup> Select strips from the pouch.
	4. Remove reagent vial caps.
	<ol> <li>Add one stirball to each new vial of Capture-R<sup>®</sup> Ready Indicator Red Cells, corQC<sup>™</sup> EXTEND Standard and DAT Positive Control Cell to be used. Gently agitate each vial to resuspend the red blood cells.</li> </ol>
	<ol> <li>Load reagents, microplates, and blood samples (and/or donor samples) onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	7. Assign the DAT assay to the blood samples (and/or donor samples), either manually or following the upload worklist procedure.
	8. Start the DAT assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the DAT assay, and records and interprets the blood sample (and/or donor samples) results.
	9. At the completion of the NEO Iris DAT assay, click the Results button on the main menu bar to access the blood sample (and/or donor samples) results.

Assay Button Abbreviation		Brief Synopsis of Assay Procedural Steps
Weak_D1	1.	Bring all reagents and blood samples to 18–30°C before testing.
Weak_D2	2.	Centrifuge the blood samples to separate the plasma from the red blood cells. Remove the caps from the blood sample tubes.
	3.	Remove the Capture-R $\ensuremath{^{ extsf{B}}}$ Select microplate frame and the desired number of Capture-R $\ensuremath{^{ extsf{B}}}$ Select strips from the pouch.
	4.	Remove reagent vial caps.
	5.	Add one stirball to each new vial of Capture-R® Ready Indicator Red Cells and corQC <sup>™</sup> EXTEND Standard to be used. Gently agitate each vial to resuspend the red blood cells.
	6.	Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.
	7.	Assign the Weak_D1 or the Weak_D2 assay to the blood samples, either manually or following the upload worklist procedure. The NEO Iris software can automatically assign the Weak_D1 or the Weak_D2 assay to the necessary blood samples if the software is configured to do so. This configuration is optional.
	8.	Start the Weak_D1 or the Weak_D2 assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the Weak_D1 or the Weak_D2 assay, and records and interprets the blood sample results.
	9.	At the completion of the NEO Iris Weak_D1 or the Weak_D2 assay, click the Results button on the main menu bar to access the blood sample results.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps		
XMatch_E XMatch_S	<ol> <li>Bring all reagents, donor segment blood samples and blood samples to 18– 30°C before testing.</li> </ol>		
, , , , , , , , , , , , , , , , , , ,	2. Centrifuge the blood sample collection tubes to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge these samples.		
	3. Remove the Capture-R <sup>®</sup> Select microplate frame and the desired number of Capture-R <sup>®</sup> Select strips from the pouch.		
	4. Remove reagent vial caps.		
	<ol> <li>Add one stirball to each new vial of Capture-R<sup>®</sup> Ready Indicator Red Cells and corQC<sup>™</sup> EXTEND Standard to be used. Gently agitate each vial to resuspend the red blood cells.</li> </ol>		
	<ol> <li>Load reagents, microplates, blood samples and donor segment blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>		
	7. Assign the XMatch assay to the blood samples, either manually or following the upload worklist procedure.		
	<ol> <li>Start the XMatch assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the XMatch assay, and records and interprets the compatibility results.</li> </ol>		
	9. At the completion of the NEO Iris XMatch assay, click the Results button on the main menu bar to access the compatibility results.		
AG_Cw	1. Bring all reagents and blood samples to 18–30°C before testing.		
AG_Fya	2. Centrifuge the blood samples to separate the plasma from the red blood cells.		
AG_Fyb	Remove the caps from the blood sample tubes.		
AG_Jka	<ol> <li>Remove the Capture-R<sup>®</sup> Select microplate frame and the desired number of Capture-R<sup>®</sup> Select strips from the pouch.</li> </ol>		
AG_JKb	4. Remove reagent vial caps.		
AG_K_II	5. Add one stirball to each new vial of Capture-R $\ensuremath{\mathbb{R}}$ Ready Indicator Red Cells and		
AG_S	corQC <sup>™</sup> EXTEND Standard to be used. Gently agitate each vial to resuspend		
AG_S_li	the red blood cells.		

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps		
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.		
	7. Assign the corresponding assay (AG_#) to the blood samples, either manually or following the upload worklist procedure. The NEO Iris software can automatically assign the corresponding assay (AG_#) to the necessary blood samples if the software is configured to do so. This configuration is optional.		
	8. Start the corresponding assay (AG_#) following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the corresponding assay (AG_#) and records and interprets the blood sample results.		
	9. At the completion of the NEO Iris corresponding assay (AG_#), click the Results button on the main menu bar to access the blood sample results.		
ANTIGENS	1. Bring all reagents and blood samples to 18–30°C before testing.		
	2. Centrifuge the blood samples to separate the plasma from the red blood cells. Remove the caps from the blood sample tubes.		
	3. 3. Remove the Capture-R® Select microplate frame and the desired number of Capture-R® Select strips from the pouch.		
	4. Remove reagent vial caps.		
	5. Add one stirball to each new vial of Capture-R® Ready Indicator Red Cells and corQC <sup>™</sup> EXTEND Standard to be used. Gently agitate each vial to resuspend the red blood cells.		
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.		
	7. Assign the ANTIGENS assay to the blood samples, either manually or following the upload worklist procedure. The NEO Iris software can automatically assign the ANTIGENS to the necessary blood samples if the software is configured to do so. This configuration is optional.		
	<ol> <li>Start the ANTIGENS following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the ANTIGENS, and records and interprets the blood sample results.</li> </ol>		

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
	9. At the completion of the NEO Iris ANTIGENS assay, click the Results button on the main menu bar to access the blood sample results.
QC	1. Bring all reagents to 18–30°C before testing.
	2. Select one clean untreated microplate for use.
	3. Remove reagent vial caps.
	<ol> <li>Add one stirball to each new vial of Referencells<sup>®</sup> A1, A2, B and 0, and corQC<sup>™</sup> EXTEND Complete to be used. Gently agitate each vial to resuspend the red blood cells.</li> </ol>
	5. Load reagents and microplates onto the NEO Iris using the procedures in the Reagent QC (see Chapter 6 – Instrument Testing Operation)
	6. The NEO Iris automatically performs the QC assay, and records and interprets the QC results.
	<ol> <li>At the completion of the NEO Iris QC assay, click the Results button on the main menu bar to access the QC results.</li> </ol>

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
CMV	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the plasma/serum from the red blood cells/clot. Remove the caps from the blood sample tubes.
	<ol> <li>Remove the Capture-CMV<sup>®</sup> microplate frame and the desired number of Capture-CMV<sup>®</sup> strips from the pouch.</li> </ol>
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-CMV Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	<ol> <li>Load reagents, microplates, and blood samples onto the NEO following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	7. Assign the CMV assay to the blood samples, either manually or following the upload worklist procedure.
	<ol> <li>Start the CMV assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO automatically performs the CMV assay, and records and interprets the blood sample results.</li> </ol>
	9. At the completion of the NEO CMV assay, press the Results button on the main menu bar to access the blood sample results.

## Test Results and Interpretation

#### Introduction

The NEO Iris generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported as Negative, Positive, Equivocal, or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

#### **Possible Test Well Results**

Possible test well results for all assays include:

Indicator	Name	Description
+	Positive	The reaction value is greater that the positive cutoff value.
-	Negative	The reaction value is less than or equal to the negative cutoff value.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater that the negative cutoff value or equal to or less than the positive cutoff value.
x	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

#### **Editing Test Results**

For selected assays, the operator can edit well results that the NEO Iris identifies as Equivocal (identified as "?"). Only equivocal results can be edited. No positive or negative results reported by the NEO Iris can be edited. Refer to the table of Result Interpretations for a list of those assays with equivocal ranges that cannot be edited.

When results are edited, the NEO Iris creates an audit trail and flags results as edited on the NEO Iris Reports. The operator is only able to edit test well results. The operator cannot edit sample result interpretations.

### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Pos	ssible Test Interpretations
ABDCHECK			
ABDCHECK2			
ABORH			
ABORH2			
ABDLONG		ABO Interpretations:	A, B, AB, O, NTD/Mixed field ?, ^INV^
ABDLONG2	+, - , ?, X	Kn Interpretations:	Positive, Negative, NTD/Mixed field ?,
ABDFULL			
ABDFULL2			
Baby_BG			
Baby_BG2			
ABOD12NC		APO Interpretations	
ABOD12_I	+, - , ?, X	ABO Interpretations:	A, B, AB, O, IND / DP, "INV"
ABODBB_I	, , ,	<u>Kn interpretations</u> .	Positive, Negative, IND / DP, "INV".
AB_CTR AB_CTR2	+, - , ?, X	ABO Interpretations:	A, B, AB, O, NTD/Mixed field ?, *INV*,
		CcEe Interpretations:	CcEe, CcEE, Ccee, CCEe, CCEE, CCee, ccEe,
			ccEE, ccee, NTD/Mixed field ?,
PHENO16	+, - , ?, X	*INV*	
PHENO12		Kell Interpretations:	K+, K- NTD/Mixed field ?, *INV*
		CcEe Interpretations:	CcEe, CcEE, Ccee, CCEe, CCEE, CCee, ccEe,
PHENO	+, - , ?, X		ccEE, ccee, IND / DP, *INV*
		Kell Interpretations:	K+, K- IND / DP, *INV*

Assay	Possible Well Results	Possible Test Interpretations		
A_subg	+, - , ?, X	A1, A2, A_int, *INV*, NTD		
AG_K Kell	+, - , ?, X	K+, K-, NTD/Mixed field ?, *INV*		
AG_M	+, - , ?, X	M+, M-, NTD/Mixed field ?, *INV*		
Cw_lgM	+, - , ?, X	Cw+, Cw-, NTD/Mixed field ?, *INV*		
AG_H	+, - , ?, X	Positive, negative, NTD, *INV*		
Pool_Cell, 2_Cell	+, -, ?*, X	Positive, Negative, No_Int, *INV*		
3_Cell Ab_ID ExtendDP, ExtendDN	+, -, ?*, X	Positive, Negative, No_Int, *INV*, Ctrl Fail Note: Ctrl Fail indicates that the per sample positive control well (4 <sup>th</sup> well in 3_Cell and 15 <sup>th</sup> well in Ab_ID, ExtendDN, and ExtendDP) has failed. This control well must react at least 3+ in order for the sample test result to be valid. If the control well is less than 3+, then the interpretation given will be Ctrl Fail.		
QC3_Cell	+,-,?,X	Qualified, FailedNote:Note:The positive control well (4th well for each set of 4 wells per test) must react at 1+ for both the weak positive control and the negative control so that the control results can be categorized as valid. If either of these 4th wells is less than 1+, then the interpretation for the QC result will be Failed. The		

Assay	Possible Well Results	Possible Test Interpretations	
		test wells for the weak positive control (1 <sup>st</sup> , 2 <sup>nd</sup> and 3 <sup>rd</sup> wells of 4) must react at 1+ so that the control result can be categorized as valid. If one or more of the 1 <sup>st</sup> , 2 <sup>nd</sup> or 3 <sup>rd</sup> wells, for the weak positive control, is less than 1+, then the interpretation for the QC result will be <b>Failed</b> .	
Weak_D1			
Weak_D2	+, -, ?, X	Positive, Negative, *INV*, NTD	
XMatch_E		Sample Interpretations: IgG Comp (Check ABO Comp),	
XMatch_S	+, - , ?*,X	No Int, *INV*	
		<u>Control Interpretations</u> : Positive, Negative, No_Int, *INV*	
DAT	+, - , ?*,X	Positive, Negative, No_Int, *INV*	
AG_Cw		Cw+, Cw-, NTD, *INV*	
AG_Fya		Fya+, Fya-, NTD, *INV*	
AG_Fyb		Fyb+, Fyb-, NTD, *INV*	
AG_Jka		Jka+, Jka, NTD, *INV*	
AG_Jkb	+,-,?*,X	Jkb+, Jkb-, NTD, *INV*	
AG_k_I	+,-,?*,X	k+, k-, NTD, *INV*	
AG_S	+,-,?*,X	S+, S-, NTD, *INV*	
AG_s_li	+,-,?*,X	s+, s-, NTD, *INV*	
ANTIGENS	+,-,?*,X	Fy(a) interpretation: Fya+, Fya-, NTD, *INV* Fy(b) interpretation: Fyb+, Fyb-, NTD, *INV* Jk(a) interpretation: Jka+, Jka, NTD, *INV* Jk(b) interpretation: Jkb+, Jkb-, NTD, *INV* S interpretation : S+, S-, NTD, *INV* s interpretation: s+, s-, NTD, *INV*	
CMV	+, - , X	Sample Interpretations: Positive, Negative, *INV* Control Interpretations: Positive, Negative, *INV*, No_Int	



**Note:** Assays that include the use of two different Anti-D reagents and yield a positive result with one Anti-D reagent, but a negative result with the other Anti-D reagent for a given sample, will generate an NTD result related to the Rh (D) testing. This can be due to, but is not limited to, reduced D antigen epitope representation on the red blood cell membrane, or the presence of a D variant. Resolution of this discrepancy can include, but is not limited to, Weak D testing and off-line testing with a battery of assorted Anti-D reagents to investigate the presence of D antigen variations, as dictated by your local on-site policies and procedures.

### Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
*	Equivocal results that may be edited. Run controls for these assays cannot be edited.
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test Interpretation	Description	
Ctrl Fail	Control Failed	
No_Int	No Interpretation	
INV	Invalid	
NTD	No Type Determined	
NTD/Mixed field ?		
IND / DP No Type Determined or possible Mixed Field		

# Attachment II: NEO Iris Operator Manual

### In This Attachment

This attachment describes the interface specifications and message structure to be used in the transfer of data from the NEO Iris to an LIS.

ATTACHMENT II: NEO IRIS OPERATOR MANUAL	II-1
Copyrights and Disclaimers	II-2
Scope	11-4
Essential Information for Communication	II-5
Result Message Structure	II-10
Host Query Message Structure	II-24
Order Message Structure	II-27
Message Examples	II-30

## Copyrights and Disclaimers

©2016, Immucor, Inc. NEO Iris™ is a trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment 2 for the NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (<u>doccontrol-de@immucor.com</u>).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## Scope

This document describes the structure to be used in the transfer of data between the Immucor NEO Iris instrument to a Laboratory Information System (LIS).

To ensure industry standardization the popular ASTM protocol is used in formatting the data structure as well as data communication transfer. This is based on the American Society for Testing and Materials designations: **E1381-02 and E1394 – 97**. These standards are more recently recognized as LIS1-A and LIS2-A, respectively.

## Essential Information for Communication

Connectivity for data transfer is accomplished with the industry standard LIS1-A (*Standard Specification for Low-Level Protocol to Transfer Messages Between Clinical Laboratory Instruments and Computer Systems* ~ Formerly ASTM E1381). Connectivity can be configured to use either Serial or Socket communication.

### **Serial Communication**

For serial communication, connectivity will make use of an available 9-pin male serial port on the NEO Iris computer. The default connection properties for serial communication are:

Baud Rate	9600
Data Bits	8
Stop Bits	1
Parity	None

The contact pin assignments for serial communication are as follows:

Direction				
Contact No.	EIA Circuit	Description	Instrument	LIS
1		Shield		No Connection
2	ВА	Receive Data	Input	Output
3	BB	Transmit Data	Output	Input
5	AB	Signal Ground		

### **Socket Communication**

For Socket communication, the traffic will traverse over a TCP/IP network, where the NEO Iris instrument will behave as the client in the connection relationship. It is expected that the LIS behave as the Server. The LIS Host/Server IP address and port number must be known prior to the NEO Iris configuration and the host must be ready (in listening mode) to establish a viable socket connection.

### Important Note Regarding Network Connectivity

When connecting to a TCP/IP network, as with a socket connection, it is important that the NEO Iris instrument be protected from malicious network threats. Risk of vulnerability is greatly minimized with the use of a firewall device. A firewall device is a hardware component, supplied by Immucor as part of the NEO Iris instrument system, that will restrict access between the NEO Iris instrument and the facility network to block unwanted use or abuse. The NEO Iris instrument will connect to the firewall device, which in turn will connect to the facility network. Serial interface connectivity does not carry through the firewall device as serial communication is not TCP/IP traffic. Please note that antivirus or additional third party software cannot be loaded on the NEO Iris computer. Doing so may impact the safety and support of the NEO Iris.

#### **General Communication Flow**

In the setup for the NEO Iris to receive order messages from the LIS, an ASTM formatted host query is generated upon loading samples on the instrument and clicking on the "Download Requests" button from the Start Run Assistant. The LIS is expected to create a response (order) to the host query message should it find matching order(s) on the LIS system. The NEO Iris is not expecting a negative or blank response if no match is found. Order messages to the NEO Iris instrument will be processed upon receipt and the assay request(s) added to the instrument Worklist.

Once results are available on the NEO Iris, the user can approve and export the results. If desired, the auto-export feature can be implemented on the NEO Iris to automatically generate the export without user intervention. This process bypasses the approval step and immediately exports results when testing is complete. Auto-export is configured per assay by Immucor Technical Support.

#### Serial and Socket Data Transfer

There are three distinct phases in transferring information between the NEO Iris instrument and the LIS. In each phase, one system directs the operation and is responsible for continuity of the communication. The three phases assure that the actions of the sender and the receiver are coordinated. The three phases are establishment, transfer, and termination. Although they are briefly described here as an overview, the LIS1-A specification should be reviewed for complete details.

#### Establishment Phase (Link Connection)

The establishment phase determines the direction of information flow and prepares the receiver to accept information. The NEO Iris will notify the LIS receiver that information is available under the following circumstances: host query and result export.

When either a host query or result message is available, the NEO Iris will initiate the establishment phase to notify the LIS that information is available. After the NEO Iris determines that the link is in a neutral state, it will transmit the <ENQ> transmission control character to the LIS. The NEO Iris will ignore all responses other than <ACK>, <NAK>, and <ENQ>. If the LIS responds with the <ACK> transmission control character, the establishment phase ends and the transfer phase begins.

When the LIS as a sender is ready to transmit order information to the NEO Iris in response to a host query, it is expected that the LIS will transmit the <ENQ> transmission control character. The NEO Iris will respond with the <ACK> transmission control character to move communication to the transfer phase.

#### **Transfer Phase**

During the transfer phase, the NEO Iris may transmit messages to or receive messages from the LIS. The transfer phase will continue until all messages are transmitted. Messages are comprised of Records which are expected to be transferred in frames, with each frame containing a maximum of 247 characters. Records longer than 240 characters are divided between two or more frames. Every record must begin a new frame. A frame is one of two types—an intermediate frame or an end frame. The frame structure is illustrated below; however, further explanation on the structure of the frame, including how the frame number and checksum are calculated, can be found in the LIS1-A specification.

h	ntermediate Frame	End Frame	
<stx>FNTe Where</stx>	ext <etb>C1C2<cr><lf></lf></cr></etb>	<stx>FNText<etx>C1C2<cr><lf></lf></cr></etx></stx>	
<stx></stx>	Start of Text transmission co	ontrol character	
FN	Single digit Frame Number 0 to 7		
Text	Data Content of Message		
<etb></etb>	End of Transmission Block transmission control character		
<etx></etx>	End of Text transmission control character		
Cl	Most significant character of check sum 0 to 9 and A to F		
C2	Least significant character of checksum 0 to 9 and A to F		
<cr></cr>	Carriage Return ASCII character		
<lf></lf>	Line Feed ASCII character		

The NEO Iris follows the rules outlined for acknowledgements to each frame. An <ACK> signifies that the last frame was received successfully and indicates readiness to receive another frame. A <NAK> signifies that the last frame was not successfully received and indicates readiness to receive the frame again. An <EOT> signifies the last frame was received successfully, but an interrupt is requested.

#### Termination Phase (Link Release)

The termination phase returns the data link to the clear or neutral state. The NEO Iris will notify the LIS that all messages have been sent by sending the <EOT> transmission control character. At this point, the NEO Iris will regard the data link to be in a neutral state. Upon receiving <EOT>, the LIS should also regard the data link to be in the neutral state.

#### Data Characters

All data will be represented as eight-bit, single-byte, coded graphic character values as defined in ISO 8859-1:1987. The eight-bit values, within the range from 0 to 127 of ISO 8859-1987 correspond to the ASCII standard character set. Values from 0 to 31 are disallowed with the exception of 7 (BEL), 9 (Horizontal Tab), 11 (Vertical Tab), and 13 (CR), where 13 is reserved as a record terminator. Values from 32 to 126 and from 128 to 254 are allowed. Values 127 and 255 are also not allowed.

Allowed characters: 7, 9, 11, 12, 13, 32-126, 128–254

Disallowed characters: 0-6, 8, 10, 14-31, 127, 255

All records are to be terminated with character 13 (CR).

The message structure and the record structure are detailed in the sections that follow.

## Result Message Structure

The following outlines the structure of result messages as they would be transmitted from the NEO Iris instrument. Result messages will contain the Header (H), Patient (P), Order (O), Result (R), and Terminator (L) records. Result messages for the IgG\_XM assay will contain the Comment (C) record to supply the identification of the donor sample used in testing.

#### **Header Record**

Bolded items are fixed data, sub-components are numbered under 'Use' column.

Field	Field Name	Use	Description of use
1	Record type ID	н	Identifies record as a header record
2	Delimiter Definition	₩^&	Definition of delimiters
3	Message Control ID	unused	
4	Access Password	unused	
5	Sender Name or ID	NEO	Fixed instrument ID
6	Sender Street Address	unused	
7	Reserved Field	unused	
8	Sender Telephone	unused	
9	Characteristic of Sender	unused	
10	Receiver ID	LIS	Fixed Receiver ID
11	Special Instructions	either Manual Edit or (null)	Test edited flag Manual Edit = test data has been edited Can be null
12	Processing ID	unused	
13	Version No.	nn	Reflects version level of the specification

14	Date and time of message	YYYYMMDDHHMMSS	Current date and time of message
----	--------------------------	----------------	----------------------------------

#### Example:

H|\^&||NEO||||LIS||LIS2-A220100219083911
#### **Patient Record**

The Patient (P) Record is sent as a placeholder. The P Record will be blank.

Field	Field Name	Use	Description of use
1	Record Type ID	Р	Identifies Record as a Patient Information record
2	Sequence Number	1	Fixed sequence number

Example:

P|1

### Order Record

Bolded items are fixed data, sub-components are numbered under 'Use' column.

Field	Field Name	Use		Description of use	
1	Record type ID	0		Identifies record as an order record	
2	Sequence Number	1		Fixed sequence number	
3	Specimen ID	variable		Sample Number	
4	Instrument Specimen ID	unused			
5	Universal Test ID	1	Universal Test ID		
		2 Universal Test ID Name			
		3 Universal Test ID Type		^^^Assay Code	
		4 Instrument Code			
6-31	Fields 6-31 are unused				

#### Example:

0|1|002650||^^^ABORH

#### **Result Record**

Bolded items are fixed data, sub-components are numbered under 'Use' column.



Note: Please refer to Table 1 for a listing of possible values for Field 3.4 – Assay Codes.



<u>Note</u>: Please refer to Table 2 for a listing of possible values for Field 4 (Measurement Value) of the Result Record.



Note: Refer to Table 3 for well identification of the result pattern.

Field	Field Name	Use		Description of use	
1	Record type ID		R	Identifies record as a result record	
2	Sequence Number		1	Fixed sequence number	
			Universal Test ID	Reserved by standard	
2	Universal Test ID	2 Universal Test ID Name		Reserved by standard	
5			Universal Test ID Type	Reserved by standard	
			variable	Instrument Assay Code	
		1 variable		Graded Well Reaction Pattern	
4	Measurement value		variable	Interpretation	
5	Units	unused			
6	Reference Range	unused			
7	Result Abnormal Flags		unused		
8	Nature of Abnormality		unused		
9	Result Status	F		Fixed – denotes final result	
10	Date of Change in Instrument	unused			
11	Operator Identification	1	variable	Performed by User Name	
		2 variable		Exported by User Name	

Field	Field Name		Use	Description of use
12	Date/Time test started		unused	
13	Date/Time test completed	YY	YYMMDDHHMMSS	Date and time created
14	Instrument Identification	1 variable		Instrument Serial Number
		2	variable	Plate Name

#### Example:

R|1|^^^ABORH|--44-33^O Positive|||||F||Melissa^Brent||20100219083911| 5030090002^UA5288913

## Table 1 - Instrument Assay Codes

	Description	Instrument Codes				
Field		Blood Grouping	Antibody Screen/Panels	Antigen Phenotyping	Other	
O Record 5.4 R Record 3.4	Assay Name	ABORH ABORH_2 ReflexABO FWD_ABORH ReflexFWD Weak_D Weak_D_F Rev_ABO ABORH_AB ABO_AB_2 RfxABO_AB FWDABO_AB FWDABO_AB RfxFWD_AB Wk_D_AB Wk_D_F_AB QCTEST_AB	Pool_Cell 2_Cell 3_Cell Ab_ID ExtendDN ExtendDP QC3_Cell	Ag_CcEe Ag_C RH2 Ag_C RH2 Ag_C RH4 Ag_E RH3 Ag_e RH5 Ag_Kell QC_CcEe QC_CCE QC_C RH2 QC_C RH2 QC_E RH3 QC_E RH3 QC_e RH5 QC_Kell	DAT CMV IgG_XM	

### Table 2 - Result Record Field 4 Measurement Values

SUB-COMPONENT 1							
Assay	Possible Values for Reaction Pattern						
All Assays	-, 1, 2, 3, 4, ?, X						
SUB-COMPONE	NT 2						
	Interpretation Values						
Assay	ABO Portion	Rh Portion					
ABORH							
ABORH_2		Positive, Negative, NTD, *INV*					
ABORH_AB							
ABO_AB_2	A, d, Ad, O, NTD, ANV						
FWD_ABORH							
FWDABO_AB							
		Positive, Negative, Pending, NTD, *INV*, Offline,					
ReflexABO	A, B, AB, O, NTD, *INV*, Offline,	RPT					
ReflexFWD	ReflexABO Pending or ReflexFWD	Weak_D Pending or Weak_D_F Pending or					
RfxABO_AB	Pending or RfxABO_AB Pending or	ReflexABO Pending or ReflexEWD Pending or					
RfxFWD_AB	RfxFWD_AB Pending	RfxABO_AB Pending or					
		RfxFWD_AB Pending					
Rev_ABO	A, B, AB, O, NTD, *INV*						

SUB-COMPONENT 2				
Assay	Interpretation Values			
Weak_D				
Weak_D_F	Positive, Negative, NTD, *INV*			
Wk_D_AB				

SUB-COMPONENT 2						
Assay	Interpretation Values					
Wk_D_F_AB						
Pool_Cell						
2_Cell	Positive, Negative, No_Int, *INV*					
DAT						
3_Cell						
Ab_ID	Desitive Negative No. Int. Ctrl Epil. *INIV/*					
ExtendDN	rositive, negative, no_int, Ctil Fail, Inv					
ExtendDP						
QCTEST						
QCTEST_AB						
QC3_Cell						
QC_CcEe						
QC_CcE	Qualified Failed					
QC_C RH2						
QC_c RH4						
QC_E RH3						
QC_e RH5						
QC_Kell						
lgG_XM	IgG Comp (Check ABO Comp), Incompatible, No_Int, *INV*					
CMV	Positive, Negative, *INV*					
Ag_C RH2	C+, C-, NTD, *INV*					
Ag_c RH4	c+, c-, NTD, *INV*					
Ag_E RH3	E+, E-, NTD, *INV*					

SUB-COMPONENT 2					
Assay	Interpretation Values				
Ag_e RH5	e+, e-, NTD, *INV*				
Ag_Kell	K+, K-, NTD, *INV*				
	C c Interpretations	E e Interpretations			
Ag_CcEe	C+ c+, C+ c-, C- c+, C+ NTD, C- NTD, NTD c+, NTD c-, NTD NTD, *INV*	E+ e+, E+ e-, E- e+, E+ NTD, E- NTD, NTD e+, NTD e-, NTD NTD, *INV*			
Assay	C c Interpretations	E Interpretations			
Ag_CcE	C+ c+, C+ c-, C- c+, C+ NTD, C- NTD, NTD c+, NTD c-, NTD NTD, *INV*	E+, E-, NTD, *INV*			

### Table 3 - LIS Well Result Identification for Result Record Field 4

### **Measurement Values**

Assay	Test	phase		Assay	Test	phase
	1	Anti-A			1	Cell 1
	2	Anti-B			2	Cell 2
ABORH	3	Anti-D series 4			3	Cell 3
ABORH_2	4	Anti-D series 5			4	Cell 4
ReflexABO	5	Monoclonal Control			5	Cell 5
	6	A1-Cell			6	Cell 6
	7	B-Cell			7	Cell 7
	1	Anti-A		Ab_ID	8	Cell 8
	2	Anti-B		ExtendDN ExtendDP	9	Cell 9
	3	Anti-A,B			10	Cell 10
ABORH_AB	4	Anti-D Series 4			11	Cell 11
ABO_AB_2 RfxABO_AB	5	Anti-D Series 5			12	Cell 12
	6	A1-Cell			13	Cell 13
	7	B Cell			14	Cell 14
	8	Monoclonal Control			15	Positive Control
	1	Anti-A			16	Negative Control
FWD_ABORH	2	Anti-B			1	Anti-C
ReflexFWD	3	Anti-D series 4			2	Anti-c
	4	Monoclonal Control		Ag_CcEe	3	Anti-E
FWDABO_AB	1	Anti-A			4	Anti-e
RfxFWD_AB	2	Anti-B			5	Monoclonal Control

Assay	Test	phase		Assay	Test	p
	3	Anti-A,B			1	Γ
	4	Anti-D Series 4		Ag_C RH2	2	
	5	Monoclonal Control			1	
	1	A1-Cell		Ag_c RH4	2	
Rev_ABO	2	B-Cell			1	
Weak D	1			Ag_E RH3	2	
Weak_D_F		Monocional Control			1	
Wk_D_AB	2			Ag_e RH5	2	
Wk_D_F_AB	2	Anti-D Series 4			1	Ì
Pool_Cell	1	Cell 1		Ag_Kell	2	
	1	Cell 1			1	,
2_Cell	2	Cell 2			2	
	1	Cell 1			3	
	2	Cell 2		4	,	
S_Cell	3	Cell 3		5		
	4	Positive Control		6		
lgG_XM	1	IgG Compatibility			7	
DAT	1	DAT result			8	
CMV	1	CMV result			1	
	1	Anti-A & A1-Cell		QC_Kell	2	
QCTEST	2	Anti-A & B-Cell			1	
	3	Anti-B & B-Cell	QCTEST_AB		2	[
	4	Anti-B & A1-Cell			3	
	5	Anti-D series 4 & QC Cell			4	

Assay	Test	phase
	1	Monoclonal Control
Ag_C RH2	2	Anti-C
	1	Monoclonal Control
Ag_c RH4	2	Anti-c
	1	Monoclonal Control
Ag_e RH3	2	Anti-E
	1	Monoclonal Control
Ag_e RH5	2	Anti-e
A 17 11	1	Monoclonal Control
Ag_Kell	2	Anti-Kell
	1	Anti-C / C positive cells
	2	Anti-c / c positive cells
	3	Anti-E / E positive cells
	4	Anti-e / e positive cells
QC_CEE	5	Anti-C / C negative cells
	6	Anti-c / c negative cells
	7	Anti-E / E negative cells
	8	Anti-e / e negative cells
	1	Anti-Kell / K positive cells
QC_Kell	2	Anti-Kell / K negative cells
	1	Anti-A & A1 Cell
OCTEST AD	2	Anti-A & B Cell
QUIESI_AB	3	Anti-B & B Cell
	4	Anti-B & A1 Cell

#### Attachment II: NEO Iris Operator Manual

Anti-A,B & A1 Cell

Anti-A,B & B Cell

Anti-A,B & QC Cell

9

10

11

Assay	Test phase			Assay	Test	phase
	6	Anti-D series 4 & A1-Cell			5	Anti-D series 4 & QC Cell
	7	Anti-D series 5 & QC Cell			6	Anti-D series 4 & A1 Cell
	8	Anti-D series 5 & B-Cell			7	Anti-D series 5 & QC Cell
					8	Anti-D series 5 & B Cell

Assay	Test phase			
	1	Anti-C		
	2	Anti-c		
Ag_CCE	3	Anti-E		
	4	Monoclonal Control		
	1	Anti-C / C positive cells		
	2	Anti-C / C negative cells		
	3	Anti-c / c positive cells		
	4	Anti-c / c negative cells		
	5	Anti-E / E positive cells		
	6	Anti-E / E negative cells		
	1	Anti-C / C positive cells		
QC_C KHZ	2	Anti-C / C negative cells		
	1	Anti-c / c positive cells		
QC_C КП4	2	Anti-c / c negative cells		
	1	Anti-E / E positive cells		
QC_E KH3	2	Anti-E / E negative cells		
	1	Anti-e / e positive cells		
UC_6 KH2	2	Anti-e / e negative cells		

#### **Comment Record**

The Comment (C) Record is sent with the IgG\_XM assay. The purpose of the C Record is to contain the ID of the donor unit used in the IgG\_XM assay.

Bolded items are fixed data, sub-components are numbered under 'Use' column.

Field	Field Name	Use		Description of use
1	Record type ID	С		Identifies record as a comment record
2	Sequence Number	1		Fixed sequence number
3	Comment source	Ι		Fixed text
4	Comment Tout	1	Donor	Fixed text
4		2	variable	Donor ID used in IgG_XM assay
5	Comment Type		unused	

Example:

C|1|I|Donor^W037908150261

#### **Terminator Record**

Bolded items are fixed data.

Field	Field Name	Use	Description of use
1	Record type ID	L	Identifies record as a terminator
2	Sequence Number	1	Fixed at 1
3	Termination Code	N	Fixed at N – normal termination

Example:

 $L \mid 1 \mid N$ 

## Host Query Message Structure

The following outlines the structure of host query from the NEO Iris instrument. Host Query messages will contain Header (H), Request (Q), and Terminator (L) Records.

#### Header Record

Bolded items are fixed data.

Field	Field Name	Use	Description of use
1	Record type ID	н	Identifies record as a header record
2	Delimiter Definition	₩^&	Definition of delimiters
3	Message Control ID	unused	
4	Access Password	unused	
5	Sender Name or ID	BBX	Fixed text
6	Sender Street Address	unused	
7	Reserved Field	unused	
8	Sender Tel.	unused	
9	Characteristic of Sender	unused	
10	Receiver ID	AURORA	Fixed text
11	Special Instructions	unused	
12	Processing ID	unused	
13	Version No.	nn	Reflects version level of the specification
14	Date and time of message	YYYYMMDDHHMMSS	Date and time order generated

#### Example:

H|\^&|||BBX|||||AURORA|||LIS2-A2|20100219083911

#### **Request Record**

Bolded items are fixed data.

Field	Field Name	Use	Description of use
1	Record type ID	Q	Identifies record as a request record
2	Sequence Number	1	Query Sequence number 1 – n
3	Starting Range ID No.	aaannn	Multiple sample barcodes separated with repeat delimiter.
4	Ending Range ID No.	unused	
5	Universal test ID	ALL	Will always be requesting all tests requested for a sample.
6	Nature of request time limit	unused	
7	Beginning request date/time	unused	
8	Ending request date/time	unused	
9	Requesting Physician name	unused	
10	Requesting Physician telephone	unused	
11	User Field 1	unused	
12	User Field 2	unused	
13	Request information status code	0	Test order request

#### Example:

• Single

Q|1|W126987||ALL|||||||0

• Multiple

Q|1|0790048\\0790032\0790016\0790003||ALL|||||||0

### **Terminator Record**

Bolded items are fixed data.

Field	Field Name	Use	Description of use
1	Record type ID	L	Identifies record as a terminator
2	Sequence Number	1	Fixed at 1
3	Termination Code	N	Fixed at $N = Normal termination$

Example:

L|1|N

## Order Message Structure

The following outlines the structure of order messages to the NEO Iris instrument. Order messages shall contain Header (H), Patient (P), Order (O), and Terminator (L) Records.

#### **Header Record**

Bolded items are fixed data.

Field	Field Name	Use	Description of use
1	Record type ID	н	Identifies record as a header record
2	Delimiter Definition	₩^&	Definition of delimiters
3	Message Control ID	unused	
4	Access Password	unused	
5	Sender Name or ID	LIS	Fixed text
6	Sender Street Address	unused	
7	Reserved Field	unused	
8	Sender Tel.	unused	
9	Characteristic of Sender	unused	
10	Receiver ID	BBX	Fixed text
11	Special Instructions	unused	
12	Processing ID	unused	
13	ASTM Version No.	nn	Aurora file version number
14	Date and time of message	YYYYMMDDHHMMSS	Date and time order generated

#### **Patient Record**

At minimum, an empty Patient (P) Record should be sent to maintain the hierarchal message structure. The LIS can optionally send data in the P Record, however, data sent in the P Record is currently ignored by NEO.

Field	Field Name	Use	Description of use
1	Record Type ID	Р	Identifies Record as a Patient Information record
2	Sequence Number	1	Sequence number (1-n)

#### **Order Record**

Bolded items are fixed data, sub-components are numbered under 'Use' column.



**Note**: Please refer to Table 1 for a listing of possible values for Assay Codes.

Field	Field Name		e	Description of use
1	Record type ID		0	Identifies record as an order record
2	Sequence Number		1	Sequence number (1-n)
3	Specimen ID	1	aaannnn	Alphanumeric Sample ID
		2	aaannnn	Alphanumeric Donation ID when
				ordering crossmatch assay.
				NULL if not crossmatch assay
4	Instrument Specimen ID		unused	
5	Universal Test ID	1	Universal Test ID	Reserved by standard
		2	Universal Test ID Name	Reserved by standard
		3	Universal Test ID Type	Reserved by standard
		4	Instrument Code	Assay Code for the order
6	Priority		R	Fixed at R = Routine
7	Request Order Date		unused	
8	Specimen Collection Date		unused	
9	Collection End Time		unused	
10	Collection Volume		unused	
11	Collector ID		unused	
12	Action Code		unused	
13	Danger Code		unused	
14	Relevant Clinical Info	unused		
15	Date Specimen Received		unused	
16	Specimen Descriptor		a	S = Sample

Field	Field Name	Use	Description of use				
			C = Crossmatch Request				
17	Ordering Physician	unused					
18	Physician Tel No.	unused					
19	User Field 1	unused					
20	User Field 2	unused					
21	Laboratory Field 1	unused					
22	Laboratory Field 2	unused					
23	Report Date/Time	unused					
24	Instrument Charge	unused					
25	Instrument Section ID	unused					
26	Report Type	F	Fixed at F = Final				
	Remaining Fields are unused.						

#### **Terminator Record**

Field	Field Name	Use	Description of use
1	Record type ID	L	Identifies record as a terminator
2	Sequence Number	1	Sequence number
3	Termination Code	N	Fixed at N = Normal termination

## Message Examples

The following illustrates examples of each of the three message types, Result, Host Query, and Order messages.

#### **Result Messages:**

#### ABORH

```
H|\^&|||NEO|||||LIS|||LIS2-A2|20100224194036
P|1
O|1|R142960||^^^ABORH
```

```
R|1|^^^ABORH|--44-33^0
Positive||||F||Donna^Brent||20100216151816|5030090012^UA5645409
L|1|N
```

#### lgG\_XM

```
H|\^&|||NEO|||||LIS|||LIS2-A2|20100224194400
P|1
O|1|R38536||^^^IgG_XM
R|1|^^^IgG_XM|-^IgG Comp (Check ABO Comp)|||||F||Jimmy^Brent||20100202092230|5030090012^SC07901917
C|1|L|Donor^LS061504
L|1|N
```

### 2\_Cell

```
H|\^&|||NEO|||||LIS|||LIS2-A2|20100224194255
P|1
0|1|24531R552567842999||^^^2_Cell
R|1|^^^2_Cell|41^Positive|||||F||Stacy^Melissa||20100209103925|5030090012^X25702
688
L|1|N
```

### FWD\_ABORH

```
H|\^&|||NEO|||||LIS|||LIS2-A2|20100224194504
P|1
O|1|R112196||^^^FWD_ABORH
R|1|^^^FWD_ABORH|-44-^B
Positive|||||F||Yolanda^Donna||20100202092221|5030090001^UA4431019
L|1|N
```

## Host Query Message:

```
H|\^&|||BBX|||||AURORA|||LIS2-A2|20100219083911
Q|1|Sample01\Sample02\Barcode0815\12345||ALL||||||||0
L|1|N
```

## **Order Messages:**

#### Non-Crossmatch

#### **Multiple Orders**

#### Crossmatch

```
H|\^&|||LIS|||||BBX|||LIS2-A2|20100219083911
P|1
O|1|107216^GC18201||^^^IgG_XM|R|||||||||C|||||||F
L|1|N
```

# **Attachment III: NEO Iris Operator Manual**

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results and interpretations, and interface data used for the ABO Titration assays.

ATTACHMENT III: NEO IRIS OPERATOR MANUAL	III-1
Copyrights and Disclaimers	
How this Attachment is Organized	111-4
ABO Titration Assay Descriptions	III-6
ABO Titration Cutoffs and Reagent Components	III-12
ABO Titration Procedural Steps	III-13
ABO Titration Results and Interpretations	III-16
Interface Specification Information	III-20

## Copyrights and Disclaimers

©2017, Immucor, Inc. NEO Iris™ is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment III for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (<u>doccontrol-de@immucor.com</u>).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

#### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

#### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (III) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (AIII) in parentheses. The first three characters (NEO Iris) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

### Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

## Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to				
$\mathbf{\Lambda}$	Warning	Potentially damaging or dangerous				
		outcomes if certain critical				
		procedural steps are ignored or				
		incorrectly executed.				
1 Ti	Consult instructions for					
	use					

## ABO Titration Assay Description

#### Introduction

#### Intended Use

The assays are for the determination of titer of IgM or IgG Anti-A and Anti-B antibodies from serum or plasma on the NEO Iris<sup>™</sup> automated analyzer.

The clinical interpretation and clinical cut offs will be established by the user.

The titer screen assays are for the detection of high titer anti-A and Anti-B antibodies from serum or plasma on the NEO Iris<sup>™</sup> automated analyzer.

#### **Description of assays**

The following table lists the ABO titration assays on the NEO Iris. The titration assays for the determination of titer of IgM antibodies exist in a horizontal and in a vertical configuration.

Antibody determination	Assay Description	Assay Button Abbreviation	Assay titer range	Samples on 96-well plate	Microplates Used	Assay configuration		
lgG	low titer range assays 8-well titration	LT_IgG_A1 LT_IgG_A2 LT_IgG_B	1 – 128 (doubling dilutions)	12				
	high titer range assays 8-well titration	HT_lgG_A1 HT_lgG_A2 HT_lgG_B	16 and 64 - 4096 (doubling dilutions from 64)	12	Capture-R® Select	Vertical		
	titer screen assays 4-well titration	screen says well ration	16 - 128 (doubling dilutions)	24				

Antibody determination	Assay Description	Assay Button Abbreviation	Assay titer range	Samples on 96-well plate	Microplates Used	Assay configuration	
IgM	low titer range assays 8-well titration	LT_lgM_A1 LT_lgM_A2 LT_lgM_B	1 - 128 (doubling dilutions)	12			
	<b>titer screen</b> assays 4-well titration	Ts_lgM_A1 Ts_lgM_A2 Ts_lgM_B	16 - 128 (doubling dilutions)	24	Galileo		
	low titer range assays 8-well titration	LT_IgM_A1 LT_IgM_A2 LT_IgM_B	1 - 128 (doubling dilutions)	8	Microplate		
	<b>titer screen</b> assays 4-well titration	Ts_lgM_A1 Ts_lgM_A2 Ts_lgM_B	16 - 128 (doubling dilutions)	24		Horizontal	

## Limitations, Warnings and Notes

lcon	Description
	Limitation:
	Capture-R® Select plates: Strips removed from pouches should be used within 16
	hours.
	Capture-R® Ready Indicator Red Cells should be used no more than 24 hours after
	a stir ball has been added to the vial. Vials of reagents other than Indicator red
	cells, that have remained continuously on the NEO Iris for 72 hours (3 days) should
	be removed and replaced with fresh vials. Vials of reagents, other than Indicator
	red cells, that are removed from the NEO Iris when not in use and refrigerated can
	be used up to their expiry dates.

lcon	Description
	<u>Warning</u> : Inspect all reagents for the presence of foam before placing on the instrument when performing any assay.
	Do not vigorously agitate reagents. Shaking will produce foam in the vial that can cause the Liquid Level Detection (LLD) function of the pipetting system to erroneously aspirate foam and/or air rather than reagent. This will produce incorrect results or an error.
	<u>Warning</u> : If you are using two or more instruments, then the specific reagent vials for each instrument must be dedicated for use on that single instrument to ensure correct reagent volume tracking. If the actual reagent volume (less than the software numeric volume) is not sufficient for the number of tests scheduled, the instrument will produce invalid results and samples will need to be rescheduled for testing.
	<u>Warning</u> : If you do not add the stirballs to the Capture-R® Indicator Red Cells suspensions or to the Referencells®, the results may be invalid or incorrect. Do not touch the stirballs. You should add them directly to the cellular reagent vials using the dispenser provided. Contamination of cellular reagents can occur if the stirballs are touched. You must only add one stirball per vial of cellular reagent. Do not add more than one stirball per vial.
	<u>Warning</u> : In rare occasions a possible result for a sample tested on a titration assays is "check for inconsistent grading". The instrument will not provide a titer result. Inconsistent grading implies that one or more well grading result within the dilution series of a sample are higher than the previous well grading result despite higher dilution (refer to table "List of Cutoffs). A sample that shows the grading failure needs to be repeated.

lcon	Description
	Warning: Serum and CPDA specimens
	For low titer IgG assays, the undiluted well of the titration series occasionally gives
	an invalid result due to the red cell monolayer being destroyed if plasma drawn in
	CPDA or serum is used. The NEO Iris will report the titration as invalid due to a
	monolayer check II failure. EDTA specimens have not shown any failures caused by
	destruction of the red cell monolayer. EDTA specimens are the preferred tube type
	to use in the low titer IgG ABO titration assays. For serum and CPDA specimens
	that give invalid results, it is recommended to repeat the test. For serum and CPDA
	specimens that persistently give invalid results, the specimen needs to be drawn in
	EDTA and the test needs to be repeated.

lcon				Description	
	Warnir specim If a ne IgM as Ts_IgM assays describ	ng: A neg nen was r gative res ssays (LT_ 1_A2, Ts_I (HT_IgG_ ped in the	ative result for the not added or that t sult or a titer belov IgM_A1, LT_IgM_A2 gM_B, Ts_IgG_A1, 1 _A1, HT_IgG_A2, HT e following table:	undiluted specimen here is no antibody o v the assay's range is 2, LT_IgM_B), titer scr Ts_IgG_A2, Ts_IgG_B) <sup>-</sup> _IgG_B) the instructio	may indicate that the detected. obtained in low titer range een assays (Ts_IgM_A1, or high titer range IgG ons on how to proceed are
		Assay	Low Titer Range	High Titer Range	Titer Screen
		lgG	The assay includes a color check as a process control for specimen addition	If result is below the range of the assay or negative, confirm with corresponding low titer range assay	If result is below the range of the assay or negative, confirm with corresponding low titer range assay
	IgM	If result is negative, confirm with appropriate reverse typing method (i.g. tube test)	n/a	If result is below the range of the assay or negative, confirm with corresponding low titer range assay	
	<u>Warnir</u> times o sample	ng: To pre over an e es insteac	event contaminatio extended period of d of loading one al	n of specimens that time it is recommend iquot serveral times o	are to be tested multiple ded to aliquot and freeze on and off the instrument.

lcon	Description
	Limitation: Samples that exhibit excessive hemolysis or lipemia, or are icteric, should not be tested on the instrument. Samples that exhibit a hemolysis concentration of more than 116 mg/dl in Low titer IgM assays, 300 mg/dl in Low titer IgG assays and 583 mg/dl in high titer IgG assays must not be tested on the instrument, because they may generate erroneous results. Color check in low titer IgG range assays can detect excessive hemolysis and invalids the result. Samples that exhibit a triglyceride concentration of more than 266 mg/dl in Low titer IgM assays, 270 mg/dl in Low titer IgG assays and 4025 mg/dl in high titer IgG assays must not be tested on the instrument, because they may generate erroneous results. Icteric samples (conjugated bilirubin) are tested until a concentration of of 23.7 mg/dl in Low titer IgM assays , 22.5 mg/dl in Low titer IgG assays and 22.4 mg/dl in high titer IgG assays without showing erroneous results. Icteric samples (unconjugated bilirubin) are tested until a concentration of 18.2 mg/dl in Low titer IgG assays without showing erroneous results.
	Limitation: At least 500µl of plasma or serum needs to be present in a sample tube to ensure that the probe picks up plasma or serum.
	<u>Warning:</u> Although specimens stored on the clot have been shown to be stable over 14 days, it is advised to store specimens off the clot as this removes the potiential for hemolysis.
	<u>Limitation:</u> For the horizontal IgM range assay, the user needs to load a full plate and select all strips to start the assay regardless of the samples to be tested on the plate. The plate is not to be reloaded to the instrument. The event log alerts the user that empty or too dense wells are detected. The message refers to empty wells on the plate that are detected by the system and does not have impact on test results.

# ABO Titration Cutoffs and Reagent Components

## List of Cutoffs

Assay	Grade	Lower Limit	Upper Limit		
Abbreviation	Glade	>	< =		
LT_lgG_A1	0	0,0	35,0		
LT_IgG_A2	?	35,0	35,0		
LT_IgG_B	1+	35,0	50,0		
HT_lgG_A1	2+	50,0	72,0		
HT_lgG_A2	3+	72,0	90,0		
HT_IgG_B					
Ts_IgG_A1					
Ts_IgG_A2	4+	90,0	99,9		
Ts_IgG_B					

Assay	Grade	Lower Limit	Upper Limit		
Abbreviation		>	< =		
LT IgM A1	0	0,0	25,0		
LT_IgM_A2	?	25,0	25,0		
LT_IgM_B	1+	25,0	35,0		
Ts_IgM_A1	2+	35,0	50,0 80,0		
Ts_lgM_A2	3+	50,0			
Ts_IgM_B	4+	80,0	99,9		

### Assay Reagent Component Grid

Reagents and Microplates	LT_IgG_A1	LT_lgG_A2	LT_lgG_B	HT_lgG_A1	HT_lgG_A2	HT_lgG_B	Ts_lgG_A1	Ts_lgG_A2	Ts_lgG_B	LT_IgM_A1	LT_IgM_A2	LT_IgM_B	Ts_IgM_A1	Ts_IgM_A2	Ts_lgM_B
Capture-R® Select	х	х	х	х	х	х	х	х	Х						
UntreatedMicroplates (barcoded)										х	х	х	х	х	х
Reagent Red Blood cells (Referencells – Group A1)	x			x			x			x			x		
Reagent Red Blood cells (Referencells – Group A2)		x			x			x			x			x	
Reagent Red Blood cells (Referencells – Group B)			x			x			х			x			х
Capture-R® Ready Indicator Red Cells	х	х	х	х	х	х	х	х	Х						
Capture® LISS	х	х	х	х	х	х	х	х	Х						

## ABO Titration Procedural Steps

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
LT_lgM_A1	
LT_lgM_A2	1. Bring all reagents and blood samples to 18–30°C before testing.
LT_IgM_B	Remove the caps from the blood sample tubes. (100µl per sample is used for
Ts_lgM_A1	LT_IgM_A1, LT_IgM_A2, LT_IgM_B assay; 15µl per sample is used for Ts_IgM_A1,
Ts_IgM_A2	Ts_IgM_A2 or Ts_IgM_B assay)
Ts IaM B	3. Remove reagent vial caps.
	4. Add one stirball to each new vial of Referencells $\ensuremath{\mathbb{R}}$ A1, A2 and B to be used.
	Gently agitate each vial to resuspend the red blood cells.
	5. Load reagents, microplates, and blood samples onto the NEO Iris™ following the
	procedures in Chapter 6 – Instrument Testing Operation.
	6. Assign the assay LT_IgM_A1, LT_IgM_A2, LT_IgM_B, Ts_IgM_A1, Ts_IgM_A2 or
	Ts_IgM_B to the blood samples, either manually or following the upload worklist
	procedure.
	7. Start the LT_IgM_A1, LT_IgM_A2, LT_IgM_B, Ts_IgM_A1, Ts_IgM_A2 or Ts_IgM_B
	assay following the procedures in Chapter 6 – Instrument Testing Operation. The
	NEO Iris™ automatically performs the LT_IgM_A1, LT_IgM_A2, LT_IgM_B, Ts_IgM_A1,
	Ts_IgM_A2 or Ts_IgM_B assay, and records and interprets blood sample results.
	8. At the completion of the NEO Iris™ LT_IgM_A1, LT_IgM_A2, LT_IgM_B, Ts_IgM_A1,
	Ts_lgM_A2 or Ts_lgM_B assay, click the Results button on the main menu bar to
	access the blood
	sample results.
	9. Assay endpoint: This is determined automatically as the last result above the cut
	off, which is set at 25 for IgM titration assays.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
LT_lgG_A1	
LT_IgG_A2	1. Bring all reagents and blood samples to 18–30°C before testing.
LT_IgG_B	Remove the caps from the blood sample tubes. (100µl per sample is used for
HT_lgG_A1	LT_lgG_A1, LT_lgG_A2, LT_lgG_B assay; 15µl per sample is used for HT_lgG_A1,
HT_lgG_A2	HT_IgG_A2, HT_IgG_B, Ts_IgG_A1, Ts_IgG_A2 or Ts_IgG_B assay)
HT_IgG_B	3. Remove the Capture-R® Select microplate frame and the desired number of
Ts_lgG_A1	Capture-R® Select strips from the pouch.
Ts_lgG_A2	5. Add one stirball to each new vial of Capture-R® Ready Indicator Red Cells and
Ts_IgG_B	Referencells ${ m I}$ A1, A2 and B to be used. Gently agitate each vial to resuspend the
	red blood cells.
	6. Load reagents, microplates, and blood samples onto the NEO Iris™ following the
	procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the LT_IgG_A1, LT_IgG_A2, LT_IgG_B, HT_IgG_A1, HT_IgG_A2, HT_IgG_B,
	Ts_IgG_A1, Ts_IgG_A2 or Ts_IgG_B assay to the blood samples, either manually or
	following the upload worklist procedure. The NEO Iris™ software can automatically
	assign the LT_IgG_A1, LT_IgG_A2, LT_IgG_B, HT_IgG_A1, HT_IgG_A2, HT_IgG_B,
	Ts_IgG_A1, Ts_IgG_A2 or Ts_IgG_B assay to the necessary blood samples if the
	software
	is configured to do so. This configuration is optional.
	8. Start the LI_IgG_A1, LI_IgG_A2, LI_IgG_B, HI_IgG_A1, HI_IgG_A2, HI_IgG_B,
	Is_IgG_A1, Is_IgG_A2 or Is_IgG_B assay following the procedures in Chapter 6 –
	Instrument Testing Operation. The NEO Iris™ automatically performs the LT_IgG_A1,
	LI_IGG_A2, LI_IGG_B, HI_IGG_A1, HI_IGG_A2, HI_IGG_B, Is_IGG_A1, Is_IGG_A2 or
	Ts_IgG_B assay, and records and interprets the blood sample results.
	9. At the completion of the NEO Iris™ LT_IgG_A1, LT_IgG_A2, LT_IgG_B, HT_IgG_A1,
	H1_IgG_A2, HT_IgG_B, Ts_IgG_A1, Ts_IgG_A2 or Ts_IgG_B assay, click the Results
	button on the main menu bar to access the blood sample results.
	10. Assay endpoint: This is determined automatically as the last result above the cut
	off, which is set at 35 for IgG titration assays.
# ABO Titration Results and Interpretations

### Introduction

The NEO Iris<sup>™</sup>generates a result for each well read by the instrument and an interpretation of the results. The ABO Titration assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative or Invalid. These results are determined by comparing the well reaction value to assayspecific cutoff values. Assay cutoff values are listed in this attachment.

### **Possible Test Well Results**

Possible test well results for ABO Titration assay includes:

Indicator	Name	Description
+	Positive	The reaction value is greater than the cutoff value.
-	Negative	The reaction value is less than or equal to the cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

The ABO Titration assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

### **Editing Test Results**

For the ABO Titration assays, no instrument generated result can be edited.

## **Result Interpretations**

The result interpretation is based on the reaction pattern of test well results. Final titer interpretation is listed in the three tables below for each titration assay.

LT_lg0	5_A1	Dilution (per well)							
LT_lg(	5_A2								
LT_lg(	G_B								
LT_IgN	/I_A1								
LT_IgN	/I_A2	1	2	4	8	16	32	64	128
LT_IgM_B									
	Negative	-	-	-	-	-	-	-	-
	1	+	-	-	-	-	-	-	-
	2	+	+	-	-	-	-	-	-
	4	+	+	+	-	-	-	-	-
	8	+	+	+	+	_	_	_	-
Titer	16	+	+	+	+	+	-	_	-
	32	+	+	+	+	+	+	-	-
	64	+	+	+	+	+	+	+	-
	≥128	+	+	+	+	+	+	+	+
	Check for inconsistent grading								

HT_lgG_A1		Dilution (per well)							
HT_lgG_A2									
HT_lg	G_B								
		16	64	128	256	512	1024	2048	4096
	≥16	+	-	-	-	-	-	-	-
	64	+	+	-	-	-	-	-	-
	128	+	+	+	-	-	-	-	-
	256	+	+	+	+	-	-	-	-
	512	+	+	+	+	+	-	-	-
Titer	1024	+	+	+	+	+	+	-	-
	2048	+	+	+	+	+	+	+	-
	4096	+	+	+	+	+	+	+	+
	<16 or negative	-	-	-	-	-	-	-	-
	Check for inconsistent grading								

Ts_lg@	5_A1	Dilution (per well)		well)	
Ts_lg@	5_A2				
Ts_lg@	5_B				
Ts_lgN	/_A1				
Ts_lgN	/_A2				
Ts_lgN	И_В				
		16	32	64	128
	16	+	-	-	-
	32	+	+	-	-
Titer	64	+	+	+	-
	≥128	+	+	+	+
	<16 or negative	-	-	-	-
	Check for inconsistent grading				

# Interface Specification Information

ABO Titration assays information for the interface specification, detailed below, is supplemental information to **Attachment II for NEO Iris Operator Manual**. Refer to **Attachment II for NEO Iris Operator Manual** for a full description of the interface specification.

### **Result Record and Field 4 Measurement Values**

Field	Description	Possible Values				
R 3.4	Assay Name	LT_IgG_A1, LT_IgG_A2, LT_IgG_B, HT_IgG_A1, HT_IgG_A2, HT_IgG_B, Ts_IgG_A1, Ts_IgG_A2, Ts_IgG_B, LT_IgM_A1, LT IgM A 2, LT IgM B,Ts IgM A1, Ts IgM A2, Ts IgM B,				
R 4.1	Original graded result	-, 1, 2, 3, 4, X				
R4.2	Interpreted Result 1	Correspondin g Assay Possible Values				
		LT_IgG_A1 LT IgG A2				
		LT_IgG_B	1, 2, 4, 8, 16, 32, 64, ≥128, Check for			
		LT_lgM_A1	inconsistent grading, *INV*			
		LT_IgM_A2				
		LT_IgM_B				
		HT_lgG_A1	≥16, 64, 128, 256, 512, 1024, 2048,			
		HT_lgG_A2	4096, <16 or negative, Check for			
		HT_IgG_B	inconsistent grading, *INV*			
		Ts_lgG_A1				
		Ts_lgG_A2				
		Ts_IgG_B 16, 32, 64, ≥128, <16 or n				
		Ts_IgM_A1 Check for inconsistent gradi				
		Ts_lgM_A2				
		Ts_lgM_B				

# LIS Well Result Identification for Field 4 Measurement Values (R 4.1)

Assay	Те	st phase
LT_IgG_A1	1	Dilution 1:1
LT_IgG_A2	2	Dilution 1:2
LT_IgG_B	3	Dilution 1:4
LT_IgM_A1	4	Dilution 1:8
LT_IgM_A2	5	Dilution 1:16
LT_IgM_B	6	Dilution 1:32
	7	Dilution 1:64
	8	Dilution 1:128

Assay	Те	est phase
HT_lgG_A1	1	Dilution 1:16
HT_lgG_A2	2	Dilution 1:64
HT_lgG_B	3	Dilution 1:128
	4	Dilution 1:256
	5	Dilution 1:512
	6	Dilution 1:1024
	7	Dilution 1:2048
	8	Dilution 1:4096

Assay	Test phase				
Ts_IgG_A1	1	Dilution 1:16			
Ts_IgG_A2	2	Dilution 1:32			
Ts_lgG_B	3	Dilution 1:64			
Ts_IgM_A1	-				
Ts_IgM_A2	4	Dilution 1:128			
Ts_lgM_B					

# **Attachment IV: NEO Iris Operator Manual**

# In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results, interpretations and interface data used for the Monoclonal Jka/Jkb assay.

ATTACHMENT IV: NEO IRIS OPERATOR MANUAL	IV-1
Copyrights and Disclaimers	IV-2
How this Attachment is Organized	IV-4
Monoclonal Jka/Jkb Assay Descriptions	IV-6
Monoclonal Jka/Jkb Assay Cutoffs and Reagent Components.	IV-9
Monoclonal Jka/Jkb Assay Procedural Steps	IV-11
Monoclonal Jka/Jkb Assay Results and Interpretations	IV-122
Monoclonal Jka/Jkb Interface Specification Information	IV-14

# Copyrights and Disclaimers

©2019, Immucor, Inc. NEO Iris<sup>®</sup> is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment IV for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (<u>doccontrol-</u> <u>de@immucor.com</u>).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

# How this Attachment is Organized

### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (IV) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (AIV) in parentheses. The first set of characters (NEO Iris\_EU) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

### Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

# Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
$\mathbf{\Lambda}$	Warning	Potentially damaging or dangerous
		outcomes if certain critical
		procedural steps are ignored or
		incorrectly executed.
1 Ti	Consult instructions for	
	use	

# Monoclonal Jka/Jkb Assay Description

### Introduction

#### Intended Use

Anti-Jk<sup>a</sup> (Monoclonal) Gamma-clone<sup>®</sup> and Anti-Jk<sup>b</sup> (Monoclonal) Gamma-clone<sup>®</sup> blood grouping reagents are intended for the detection of Jka (JK1) and Jkb (JK2) antigens, respectively, on red blood cells by direct hemagglutination on the NEO Iris<sup>®</sup> automated analyzer.

#### List of assays

The following table lists the monoclonal Jka/Jkb assay on the NEO Iris. The assay exists in a horizontal and in a vertical configuration.

Due to a transition period two (2) assay versions will be available for the NEO Iris for a certain time.

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
Monoclonal Jka/Jkb	AG_Jkab_m	<ol> <li>Monoclonal Control</li> <li>Anti-Jk<sup>a</sup> (Monoclonal) Gamma- clone<sup>®</sup></li> <li>Anti-Jk<sup>b</sup> (Monoclonal) Gamma- clone<sup>®</sup></li> </ol>	Untreated microplates
	Blood grouping	Vertical/8 wells per strip	
Monoclonal Jka/Jkb	AG_Jkab_m	<ol> <li>Monoclonal Control</li> <li>Anti-Jk<sup>a</sup> (Monoclonal) Gamma- clone<sup>®</sup></li> <li>Anti-Jk<sup>b</sup> (Monoclonal) Gamma- clone<sup>®</sup></li> </ol>	Untreated microplates

#### Version 1.00:

#### Version 2.00i:

Assay Description	Assay short name		Used reagents	Microplates Used
	Blood grouping		Horizontal/12 wells per strip	
Monoclonal Jka/Jkb	AG_Jkab_m	<ol> <li>1.</li> <li>2.</li> <li>3.</li> <li>4.</li> <li>5.</li> </ol>	Monoclonal Control Anti-Jk <sup>a</sup> (Monoclonal) Gamma- clone <sup>®</sup> Anti-Jk <sup>b</sup> (Monoclonal) Gamma- clone <sup>®</sup> GammaZyme-B <sup>TM</sup> Specimen Diluent	Untreated microplates
	Blood grouping		Vertical/8 wells per strip	
Monoclonal Jka/Jkb	AG_Jkab_m	<ol> <li>1.</li> <li>2.</li> <li>3.</li> <li>4.</li> <li>5.</li> </ol>	Monoclonal Control Anti-Jk <sup>a</sup> (Monoclonal) Gamma- clone <sup>®</sup> Anti-Jk <sup>b</sup> (Monoclonal) Gamma- clone <sup>®</sup> GammaZyme-B <sup>TM</sup> Specimen Diluent	Untreated microplates

# Limitations, Warnings and Notes

lcon	Description
	Limitation: The purpose of the Monoclonal Control for AG_Jkab_m assay is to serve as a sample control for Immucor low protein blood grouping reagents. It is expected to indicate those sample related conditions that could lead to spontaneous agglutination with low protein reagents and, therefore, a false positive interpretation for the test. When the Monoclonal Control well yields a positive result, then the instrument would not report the blood type.
	Limitation: Vials of reagents, that have remained continuously on the NEO Iris for 72 hours (3 days) should be removed and replaced with fresh vials. Vials of reagents that are removed from the NEO Iris when not in use and refrigerated can be used up to their expiration dates.
	Limitation: The NEO Iris cannot reliably detect hemagglutination reactions that are graded as 1+ or less in test tube methodology. The NEO Iris does not generate an interpretation of mixed-field. Such a mixed-field reaction will be interpreted as positive, negative, or equivocal.
	Limitation: The grading of reactions on the NEO Iris must only be regarded as an approximation when compared to off-line visual grading by laboratory technical staff.

# Monoclonal Jka/Jkb Assay Cutoffs and Reagent Components

# List of Cutoffs

Version 1.00:

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		0	0	20
		?	20	23
	Monoclonal control	1+	23	35
		2+	35	50
		3+	50	75
		4+	75	99
		0	0	20
	Anti-Jk <sup>a</sup> (Monoclonal) Gamma-clone®	?	20	23
AC IIIah m		1+	23	35
		2+	35	50
		3+	50	75
		4+	75	99
		0	0	20
		?	20	23
	Anti-Jk <sup>b</sup> (Monoclonal) Gamma-clone®	1+	23	35
		2+	35	50
		3+	50	75
		4+	75	99

Version 2.00i:

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	< =
		0	0	25
		?	25	28
		1+	28	35
	Monoclonal	2+	35	50
		3+	50	80
		4+	80	99
		0	0	25
	Anti-Jk <sup>a</sup>	?	25	28
		1+	28	35
АĞ_ЈКар_Ш	Gamma-clone <sup>®</sup>	2+	35	50
		3+	50	80
		4+	80	99
		0	0	25
		?	25	28
	Anti-Jk <sup>b</sup> (Monoclonal) Gamma-clone®	1+	28	35
		2+	35	50
		3+	50	80
		4+	80	99

# Monoclonal Jka/Jkb Assay Procedural Steps

## Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button		
Abbreviation		Brief Synopsis of Assay Procedural Steps
AG Jkab m		
	1.	Bring all reagents and blood samples to 18–30°C before testing.
	2.	Centrifuge the blood samples to separate the plasma from the red blood
		cells and then remove the caps from those tubes. Process the donor
		segment blood samples, but do not centrifuge those samples.
	3.	Remove reagent vial caps.
	4.	Load reagents, microplates, and blood samples, onto the NEO Iris following
		the procedures in Chapter 6 – Instrument Testing Operation.
	5.	Assign the AG_Jkab_m assay to the blood samples either manually or
		following the upload worklist procedure.
	6. Start the AG_Jkab_m assay following the procedures in Chapter 6 –	
		Instrument Testing Operation. The NEO Iris automatically performs the
		AG_Jkab_m assay, and records and interprets blood sample results.
	7.	At the completion of the NEO Iris AG_Jkab_m assay, click the Results button
		on the main menu bar to access the blood samples results.

# Monoclonal Jka/Jkb Assay Results and Interpretations

#### Introduction

The NEO Iris<sup>®</sup>generates a result for each well read by the instrument and an interpretation of the results. The Monoclonal Jka/Jkb Assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative, Invalid or Equivocal. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

### **Possible Test Well Results**

Possible test well results for Monoclonal Jka/Jkb Assay includes:

Indicator	Name	Description	
+	Positive	The reaction value can be considered as positive.	
-	Negative	The reaction value can be considered as negative.	
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.	
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).	

## **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations
AG_Jkab_	+(1, 2, 3, 4), -, ? <sup>1</sup> , X	Kidd interpretations: Jka+, Jka-, Jkb+,
m		Jkb-, NTD, Ctrl failure <sup>2</sup> , *INV*

<sup>1</sup> Equivocal results that may be edited.

<sup>2</sup> Test Interpretation only applicable for Assay Version 2.00i

### Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test Interpretation	Description
INV	Invalid
NTD	No Type Determined
Ctrl failure*	Control failure

\* Test Interpretation only applicable for Assay Version 2.00i

# Interface Specification Information

The Monoclonal Jka/Jkb assay information for the interface specification, detailed below, is supplemental information to **Attachment II for the NEO Iris Operator Manual (EU).** Refer to **Attachment II for the NEO Iris Operator Manual (EU)** for a full description of the interface specification.

### **Result Record Field 4 Measurement Values**

SUB-COMPONENT 1				
Assay	Possible Values for Reaction Pattern			
All Assays	-, 1, 2, 3, 4, ?, X			
SUB-COMPONENT 2				
Assay	Interpretation Values			
AC Ikah m	Jka Interpretation	Jkb Interpretation		
	Jka+, Jka-, NTD, Ctrl failure <sup>1</sup> , *INV*	Jkb+, Jkb-, NTD, *INV*		
AC likeb m <sup>2</sup>	Jka+ Jkb+, Jka- Jkb+, Jka+ Jkb-, Jka- Jkb-, NTD, Ctrl failure <sup>1</sup> ,			
	*INI\/*			

<sup>1</sup> Interpretation value only applicable for Assay Version 2.00i

<sup>2</sup> Applicable for alternative Aurora file with one Element

### LIS Well Result Identification for Result Field 4 Measurement Values

Assay	Test phase	
	1	Anti-Jk <sup>a</sup> (Monoclonal) Gamma-clone®
AG_Jkab_m	2	Anti-Jk <sup>b</sup> (Monoclonal) Gamma-clone®
	3	Monoclonal Control

# **Attachment V: NEO Iris Operator Manual**

# In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results and interpretations, and interface data used for the modified Antigens Assays ANTIGENS4 and ANTIGENS5.

Copyrights and DisclaimersV-2 How this Attachment is OrganizedV-4 Modified Antigens Assays DescriptionsV-6 Modified Antigens Assays Cutoffs and Reagent ComponentsV-8 Modified Antigens Assays Procedural StepsV-9 Modified Antigens Assays Results and InterpretationsV-11	ŀ	ATTACHMENT V: NEO IRIS OPERATOR MANUAL	.V-1
How this Attachment is OrganizedV-4 Modified Antigens Assays DescriptionsV-6 Modified Antigens Assays Cutoffs and Reagent ComponentsV-8 Modified Antigens Assays Procedural StepsV-9 Modified Antigens Assays Results and InterpretationsV-11		Copyrights and Disclaimers	.V-2
Modified Antigens Assays DescriptionsV-6 Modified Antigens Assays Cutoffs and Reagent ComponentsV-8 Modified Antigens Assays Procedural StepsV-9 Modified Antigens Assays Results and InterpretationsV-11		How this Attachment is Organized	.V-4
Modified Antigens Assays Cutoffs and Reagent ComponentsV-8 Modified Antigens Assays Procedural StepsV-9 Modified Antigens Assays Results and InterpretationsV-11		Modified Antigens Assays Descriptions	.V-6
Modified Antigens Assays Procedural StepsV-9 Modified Antigens Assays Results and InterpretationsV-11		Modified Antigens Assays Cutoffs and Reagent Components	.V-8
Modified Antigens Assays Results and Interpretations V-11		Modified Antigens Assays Procedural Steps	.V-9
		Modified Antigens Assays Results and Interpretations	/-11

# Copyrights and Disclaimers

©2017, Immucor, Inc. NEO Iris™ is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment V for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (<u>doccontrol-de@immucor.com</u>).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

# How this Attachment is Organized

### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (V) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (AV) in parentheses. The first set of characters (NEO Iris\_EU) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

## Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

# Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
$\mathbf{\Lambda}$	Warning	Potentially damaging or dangerous
		outcomes if certain critical
		procedural steps are ignored or
		incorrectly executed.
1 îi	Consult instructions for	
لطما	use	

# Modified Antigens Assays Description

#### Introduction

#### Intended Use

Within the original Antigens assay (ANTIGENS) Immucors blood grouping reagents Anti-Fy(a) micro, Anti-Fy(b) micro, Anti-Jk(a) micro and Anti-Jk(b) micro, Anti-S micro, Anti-s micro are intended for the detection of the Fya (FY1), Fyb (FY2), Jka (JK1) and Jkb (JK2) antigens, S (MNS 3), s (MNS 4) antigens on red blood cells by the Capture R method.

The two modified Antigens assays (ANTIGENS4 and ANTIGENS5) are intended for the detection of following rare antigens on human erythrocytes:

**ANTIGENS4:** Fya (FY1), Fyb (FY2), S (MNS 3), s (MNS 4) antigens on red blood cells by the Capture R method.

**ANTIGENS5:** Fya (FY1), Fyb (FY2), Jka (JK1), S (MNS 3), s (MNS 4) antigens on red blood cells by the Capture R method.

#### List of assays

The following table lists the modified Antigens assays on the NEO Iris.

Assay Description	Assay short name	Used reagents	Microplates Used
	Capture Select -		
	Rare Antigens		
	ANTIGENS4	Anti-Fy(a) micro	
Capture-R Antigens		Anti-Fy(b) micro	
		Empty	
		Empty	Capture D Calact Distan
		Anti-S micro	Capture-R Select Plates
		Anti-s micro	
		Empty	
		Negative Control micro	

Assay Description	Assay short name	Used reagents	Microplates Used
Capture-R Antigens		Anti-Fy(a) micro	
		Anti-Fy(b) micro	
		Anti-Jk(a) micro	
		Empty	Conturo D. Coloct Distor
	AINTIGENSS	Anti-S micro	Captule-R Select Plates
		Anti-s micro	
		Empty	
		Negative Control micro	

# Limitations, Warnings and Notes

lcon	Description
	Limitation:
<b></b>	Capture-R® Ready Indicator Red Cells should be used no more than 24 hours after astir ball has been added to the vial. Vials of reagents, other than Indicator Red Cells, that have
	remained continuously on the NEO Iris for 72 hours (3 days) should be removed and replaced with fresh vials. Vials of reagents, other than Indicator Red Cells, that are removed from the NEO Iris when not in use and refrigerated can be used up to their expiration dates.
	Limitation: Capture-R® Select plates: Strips removed from pouches should be used within 16 hours.
	Limitation: The grading of reactions on the NEO Iris must only be regarded as an approximation when compared to off-line visual grading by laboratory technical staff.

# Modified Antigens Assays Cutoffs and Reagent Components

# List of Cutoffs

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		0	0	20
	Anti-Fy(a) micro	?	20	40
ANTIGENS4	Anti-Fy(b) micro Anti-S micro Anti-s micro	1+	40	50
		2+	50	72
	Negative Control micro	3+	72	90
		4+	90	100

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		0	0	20
	Anti-Fy(a) micro Anti-Fy(b) micro Anti-Jk(a) micro	?	20	40
ANTIGENS5		1+	40	50
	Anti-S micro	2+	50	72
	Anti-s micro	3+	72	90
	Negative Control micro	4+	90	100

# Modified Antigens Assays Procedural Steps

# Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
ANTIGENS4	. Bring all reagents and blood samples to 18–30°C before testing.
ANTIGENS5	<ol> <li>Centrifuge the blood samples to separate the plasma from the red blood cells.</li> <li>Remove the caps from the blood sample tubes.</li> </ol>
	8. Remove the Capture-R <sup>®</sup> Select microplate frame and the desired number of Capture-R <sup>®</sup> Select strips from the pouch.
	. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-R® Ready Indicator Red Cells to be used. Gently agitate each vial to resuspend the red blood cells.
	<ol> <li>Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	7. Assign the assay ANTIGENS4 / ANTIGENS5 to the blood samples, either manually or following the upload worklist procedure. The NEO Iris software can automatically assign the assay ANTIGENS4 / ANTIGENS5 to the necessary blood samples if the software is configured to do so. This configuration is optional.
	8. Start the assay ANTIGENS4 / ANTIGENS5 following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the assay ANTIGENS4 / ANTIGENS5, and records and interprets the blood smaple results.
	<ol> <li>At the completion of the NEO Iris assay ANTIGENS4 / ANTIGENS5, click the Results button on the main menu bar to access the blood sample results.</li> </ol>

# Modified Antigens Assays Results and Interpretations

### Introduction

The NEO Iris generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported as Negative, Positive, Equivocal, or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

### **Possible Test Well Results**

Possible test well results for the modified Antigens assays ANTIGENS4 and ANTIGENS5 include:

Indicator	Name	Description
+	Positive	The reaction value is greater than the positive cutoff value.
-	Negative	The reaction value is less than or equal to the negative cutoff value.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

## **Editing Test Results**

For selected assays, the operator can edit well results that the NEO Iris identifies as Equivocal (identified as "?"). Only equivocal results can be edited. No positive or negative results reported by the NEO Iris can be edited. Refer to the table of Result Interpretations for a list of those assays with equivocal ranges that cannot be edited.

When results are edited, the NEO Iris creates an audit trail and flags results as edited on the NEO Iris Reports. The operator is only able to edit test well results. The operator cannot edit sample result interpretations.

## **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations
ANTIGENS4	+,-,?*,X	Fy(a) interpretation: Fya+, Fya-, NTD, *INV*
		Fy(b) interpretation: Fyb+, Fyb-, NTD, *INV*
		S interpretation : S+, S-, NTD, *INV*
		s interpretation: s+, s-, NTD, *INV*
ANTIGENS5	+,-,?*,X	Fy(a) interpretation: Fya+, Fya-, NTD, *INV*
		Fy(b) interpretation: Fyb+, Fyb-, NTD, *INV*
		Jk(a) interpretation: Jka+, Jka-, NTD, *INV*
		S interpretation : S+, S-, NTD, *INV*
		s interpretation: s+, s-, NTD, *INV*

\*Equivocal results that may be edited.

# Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test Interpretation	Description
INV	Invalid
NTD	No Type Determined

# **Attachment VI: NEO Iris Operator Manual**

# In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results and interpretations, and interface data used for assays to provide basic blood type, Rh status and/or Kell and/or CDE.

ATTACHMENT VI: NEO IRIS OPERATOR MANUAL	VI-1
Copyrights and Disclaimers	VI-2
How this Attachment is Organized	VI-4
Assay Descriptions	VI-6
Cutoffs and Reagent Components	VI-12
Procedural Steps	VI-54
Results and Interpretations	VI-56

# Copyrights and Disclaimers

©2018, Immucor. Inc. NEO Iris™ is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment VI for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (<u>doccontrol-</u> <u>de@immucor.com</u>).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.
# How this Attachment is Organized

## In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

## **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (VI) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (AVI) in parentheses. The first three characters (NEO Iris) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

## Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

# Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to	

	Warning	Potentially damaging or dangerous
		outcomes if certain critical
		procedural steps are ignored or
		incorrectly executed.
1 îi	Consult instructions for	
	use	

# Assay Description

### Introduction

### Intended Use

The described assays are designed to provide basic blood type, Rh status and/or Kell and/or CDE which can include both forward testing (of red blood cell) and reverse testing (of plasma/serum). Assays are based on direct agglutination technology for the NEO Iris<sup>™</sup> automated analyzer.

### List of assays

The following table lists the described assays on the NEO Iris. The assays exist in a horizontal and/or vertical configuration.

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
Forward and Reverse ABO Blood Grouping	ABDLONG	<ol> <li>NOVACLONE Anti-A</li> <li>NOVACLONE Anti-B</li> <li>NOVACLONE Anti-AB</li> <li>A1-Cell</li> <li>A2-Cell</li> <li>B-Cell</li> <li>O-Cell</li> <li>Autocontrol</li> <li>NOVACLONE Anti-D</li> <li>immuClone Anti-D rapid</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates

### Assay Description

Assay Description	Assay short name	Used reagents	Microplates Used	
		1. immuClone Anti-A		
		2. immuClone Anti-B		
		3. immuClone Anti-AB		
		4. A1-Cell		
Forward and		5. A2-Cell		
Reverse ABO Blood	ABDLONG2	6. B-Cell	Untreated microplates	
Grouping		7. O-Cell		
		8. Autocontrol		
		9. NOVACLONE Anti-D		
		10. immuClone Anti-D rapid		
		11. immuClone Rh-Hr Control		
		1. NOVACLONE Anti-A		
		2. NOVACLONE Anti-B		
		3. NOVACLONE Anti-AB		
		4. A1-Cell		
Forward and		5. A2-Cell		
Reverse ABO Blood		6. B-Cell	Untreated microplates	
Grouping incl.	ADULUNGS	7. O-Cell	ontreated microplates	
CDE		8. Autocontrol		
		9. NOVACLONE Anti-D		
		10. immuClone Anti-D rapid		
		11. immuClone Anti-CDE		
		12. immuClone Rh-Hr Control		

Assay Description	Assay short name	Used reagents	Microplates Used	
		1. immuClone Anti-A		
		2. immuClone Anti-B		
		3. immuClone Anti-AB		
		4. A1-Cell		
Forward and		5. A2-Cell		
Reverse ABO		6. B-Cell	Untropted microplates	
Grouping incl.	ABDLONG4	7. O-Cell	Untreated microplates	
CDE		8. Autocontrol		
		9. NOVACLONE Anti-D		
		10. immuClone Anti-D rapid		
		11. immuClone Anti-CDE		
		12. immuClone Rh-Hr Control		
		1. NOVACLONE Anti-A		
		2. NOVACLONE Anti-B		
		3. NOVACLONE Anti-AB		
		4. A1-Cell		
Forward and		5. A2-Cell		
Reverse ABO		6. B-Cell	Untropted microplates	
Grouping incl.	ABDLONGK	7. O-Cell	Untreated microplates	
Kell		8. Autocontrol		
		9. NOVACLONE Anti-D		
		10. immuClone Anti-D rapid		
		11. Automated immuClone Anti-Kell		
		12. immuClone Rh-Hr Control		

### Assay Description

Assay Description Assay short name Use		Used reagents	Microplates Used	
		1. immuClone Anti-A		
		2. immuClone Anti-B		
		3. immuClone Anti-AB		
		4. A1-Cell		
Forward and		5. A2-Cell		
Reverse ABO		6. B-Cell		
Blood Grouping incl	ABDLONG2K	7. O-Cell	Untreated microplates	
Kell		8. Autocontrol		
		9. NOVACLONE Anti-D		
		10. immuClone Anti-D rapid		
		11. Automated immuClone Anti-Kell		
		12. immuClone Rh-Hr Control		
	Blood grouping	Horizontal/6 wells per strip		
	ABD6_I	1. ImmuClone Rh-Hr Control		
5 1.150		2. immuClone Anti-A		
Forward ABO		3. immuClone Anti-B	Untreated microplates	
Grouping		4. immuClone Anti-AB		
		5. immuClone Anti-D rapid		
		6. NOVACLONE Anti-D		
		1. NOVACLONE Diluent Control		
		2. NOVACLONE Anti-A	Untreated microplates	
Forward ABO		3. NOVACLONE Anti-B		
Grouping	ABD6_N	4. NOVACLONE Anti-AB		
		5. immuClone Anti-D rapid		
		6. NOVACLONE Anti-D		

#### **Assay Description**

Assay Description		Used reagents	Microplates Used	
		1. ImmuClone Rh-Hr Control		
Forward ABO		2. immuClone Anti-A		
Blood		3. immuClone Anti-B		
Grouping incl.	ABD6CDE_I	4. immuClone Anti-AB	Untreated microplates	
CDE		5. immuClone Anti-D rapid		
		6. immuClone Anti-CDE		
		1. NOVACLONE Diluent Control		
Forward ABO		2. NOVACLONE Anti-A		
Blood		3. NOVACLONE Anti-B		
Grouping incl.	ABD6CDE_N	4. NOVACLONE Anti-AB	Untreated microplates	
CDE		5. immuClone Anti-D rapid		
		6. immuClone Anti-CDE		
	ABD6K_I	1. ImmuClone Rh-Hr Control		
Forward ABO		2. immuClone Anti-A		
Blood		3. immuClone Anti-B	Untreated microplates	
Grouping incl.		4. immuClone Anti-AB		
Kell		5. immuClone Anti-D rapid		
		6. Automated immuClone Anti-Kell		
		1. NOVACLONE Diluent Control		
Forward ABO		2. NOVACLONE Anti-A		
Blood	ABD6K_N	3. NOVACLONE Anti-B		
Grouping incl.		4. NOVACLONE Anti-AB	Untreated microplates	
Kell		5. immuClone Anti-D rapid		
		6. Automated immuClone Anti-Kell		
	Blood grouping	Horizontal/4 wells per strip		
		1. A1 cell		
Reverse ABO		2. A2 cell	Untreated microplates	
Groupina		3. B cell	Untreated microplates	
r - 2		4. O cell		
	Blood grouping	Vertical/4 wells per strip		

#### **Assay Description**

Assay Description	Assay short name	Used reagents	Microplates Used
Reverse ABO Blood Grouping	ABO_REV	1. A1 cell 2. A2 cell 3. B cell 4. O cell	Untreated microplates
	Blood grouping	Vertical/8 wells per strip	
Forward and Reverse ABO Blood Grouping	ABODFULL	<ol> <li>NOVACLONE Diluent Control</li> <li>NOVACLONE Anti-A</li> <li>NOVACLONE Anti-B</li> <li>immuClone Anti-D rapid</li> <li>NOVACLONE Anti-D</li> <li>A1-Cell</li> <li>B-Cell</li> <li>O-Cell</li> </ol>	Untreated microplates
Forward and Reverse ABO Blood Grouping	ABODFULL2	<ol> <li>immuClone Rh-Hr Control</li> <li>immuClone Anti-A</li> <li>immuClone Anti-B</li> <li>immuClone Anti-D rapid</li> <li>NOVACLONE Anti-D</li> <li>A1-Cell</li> <li>B-Cell</li> <li>O-Cell</li> </ol>	Untreated microplates

## Limitations, Warnings and Notes

All further limitations, warnings and notes that relate to the described assays are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings and in Attachment I of the NEO Iris instrument.

lcon

Description

lcon	Description
	Limitation: The Reverse ABO Blood Group assay (i.e. ABO_REV) is intended to
	screen for ABO antibodies in plasma. This provides a presumptive indication of the
	red blood cell ABO group for the individual from which the plasma was collected.
	The absence of the corresponding forward group test subjects the results of this
	assay to several potential sources of error. These include serological factors due to
	the donor's age, immunological state, the presence of an IgM allo- or
	autoantibody, transfusion or transplantation status or disease states. Sample
	related conditions such as higher-than-expected levels of lipids, bilirubin, free
	plasma hemoglobin, clots or aggregates may also be a factor. This assay cannot be
	used (alone) to determine the ABO group as part of pre-transfusion testing in a
	patient population, or to determine ABO group for blood components intended for
	transfusion.

# Cutoffs and Reagent Components

# List of Cutoffs

Assay	Postion	Crada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
	NOVACLONE Anti-A	0	0	30
		?	30	58
		1+ (not reported)	58N/A	58N/A
ABDLONG		2+ (not reported)	N/A58	N/A58
		3+ (not reported)	N/A58	N/A58
		4+	58	100
	NOVACLONE	0	0	30
	Anti-B	?	30	76

Assay	Reaction	Grade	Lower Limit	Upper Limit
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
	NOVACIONE	1+ (not reported)	N/A	N/A
	Anti-AB	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
	A1-Cell	0	0	23
		?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
		2+	35	50

Assay	Postion	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	B-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
	O-Cell	?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	Autocontrol	2+	35	50
		3+	50	76
		4+	76	100
		0	0	30
	NOVACLONE	?	30	76
	Anti-D	1+ (not reported)	N/A	N/A

Assay	Peaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		? 30	76	
		1+ (not reported)	N/A	N/A
	Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
	ImmuClone Rh-Hr Control	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay	Postion	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	< =

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		0	0	30
		?	30	58
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
		?	30	76
		1+ (not reported) N/A	N/A	
ADDLONGZ	ImmuClone Anti-B	2+ (not reported)	N/A	N/A
		3+ (not reported)	2+ (not     N/A     N/A       reported)     3+ (not     N/A     N/A       reported)     N/A     N/A	N/A
	4+	4+	76	100
		0	0	30
		?	30	76
		$ \begin{array}{c c c c c c c c c c c c c c c c c c c $	N/A	
	Anti-AB	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay	Postion	Crada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		0	0	23
		?	23	28
		1+	28	35
	AT-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	Az-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
	P. Coll	1+	28	35
	D-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	U-Cell	2+	35	50
		3+	50	76
		4+	76	100

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		0	0	23
		?	23	28
		1+	28	35
	Autocontrol	2+	35	50
		3+	50	76
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	NOVACLONE Anti-D	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
	ImmuClone Anti-D rapid	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
	ImmuClone	0	0	30
	Rh-Hr Control	?	30	76

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	< =
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	Anti-A	2+ (not reported)	N/A	N/A
ABDLONG3		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
		?	30	76
	NOVACLONE Anti-B	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-AB	2+ (not reported)	N/A N N/A N 76 1 0 2	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	AT-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	A2-Cell	2+	35	50
		3+	50	76
		4+	76	100
	R-Call	0	0	23
	D-Cell	?	23	28

Assay	say	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	0         23           23         28           28         35           35         50           50         76           76         100           0         23           23         28           23         28           35         50	35
	U-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
	Autocontrol	1+	28	35
	Autocontrol	2+	35	50
		3+	50	76
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-D	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay	Peaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Gidde	>	<=
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
	ImmuClone Anti-CDE	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
	ImmuClone Rh-Hr Control	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		0	0	30
		?	30	58
	ImmuClana	1+ (not reported)	N/A	N/A
	Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
	ImmuClone Anti-B	0	0	30
		?	30	76
ABDLONG4		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		4+       58       10         0       0       3         ?       30       7         1+ (not reported)       N/A       N/         2+ (not reported)       N/A       N/         3+ (not reported)       N/A       N/         3+ (not reported)       N/A       N/         4+       76       10	N/A	
		4+	76	100
		0	0	30
		?	30	76
	ImmuClone Anti-AB	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	AT-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	A2-Cell	2+	35	50
		3+	50	76       100         0       23         23       28         28       35         35       50         50       76         76       100         0       23         23       28
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	B-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
	U-Cell	1+	28	35
		2+	35	50

Assay	Peaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
	Autocontrol	1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	30
		?	30	76
	NOVACLONE Anti-D	1+	76	76
		2+	76	76
		3+	76	76
		4+	76	100
		0	0	30
		?	30	76
	ImmuClana	1+ (not reported)	1+ (not reported) N/A	N/A
	ImmuClone Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
	ImmuClone	0	0	30
	Anti-CDE	?	30	76

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
	ImmuClone Rh-Hr Control	?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay	Reaction	Grada	Lower Limit	Upper Limit
Abbreviation		Grade	>	< =
		0	0	30
		?	30	58
ABDLONGK	NOVACLONE	1+ (not reported)	N/A	N/A
	Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		4+	58	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-B	2+ (not reported)	N/A	N/A
	3+ (not reported) N/A	N/A		
		4+	76	100 30 76 N/A N/A N/A 100 30 76 N/A N/A N/A N/A N/A 100 23 28 35
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	NOVACLONE Anti-AB	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
	A2-Cell	0	0	23

Assay	Peaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	B-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	U-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
Autocontrol		?	23	28
	1+	28	35	
	Autocontrol	2+	35	50
		3+	50	76
		4+	76	100
	NOVACLONE	0	0	30

Assay	Peaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
	Anti-D	?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	23
		?	23	50
	Automated	1+ (not reported)	N/A	N/A
	Anti-Kell	2+ (not reported)	N/A	N/A
		3+	50	80
		4+	80	100
	ImmuClone	0	0	30

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
	Rh-Hr Control	?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit
		0	0	30
		?	30	58
ABDLONG2K		1+ (not reported)N/AN/A2+ (not reported)N/AN/A3+ (not reported)N/AN/A	N/A	
	Anti-A		N/A	
			N/A	
		4+	58	100
		0	0	30
		?	30	76
	ImmuClone Anti-B	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-AB	2+ (not reported)	N/A	N/A
		3+ (not reported) N/A	N/A	
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	AT-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
	الم2-2ما	1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
	B-Cell	0	0	23

Assay	Reaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	O-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
	Autocontrol	1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	30
		?	30	76
	NOVACLONE	1+ (not reported)	N/A	N/A
	Anti-D	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A         N/A           N/A         N/A           N/A         N/A           N/A         N/A           76         100           0         23           23         50           N/A         N/A           N/A         N/A	
		4+	76	100
	Automated	0	0	23
		?	23	50
		1+ (not reported)	N/A	N/A
	Anti-Kell	2+ (not reported)	N/A	N/A
		3+	50	80
		4+	80	100.
		0	0	30
		?	30	76
	ImmuClone	1+ (not reported)	N/A	N/A
	Rh-Hr Control	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		4+	76	100

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		Giude	>	< =
	ImmuClone Rh-Hr Control	0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
	ImmuClone Anti-A	0	0	30
ABD6_I		?	30	58
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported) N/A	N/A	
		4+	58	100
		0	0	30
	ImmuClone	?	30	76
	Anti-B	1+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		2+ (not reported)	N/A	×≡ N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-AB	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
	ImmuClone Anti-D rapid	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
	NOVACLONE	?	30	76
	Anti-D	1+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	< =
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
	NOVACLONE Diluent Control	0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
ABD6_N		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
	NOVACLONE	0	0	30

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit
	Anti-B	?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-AB	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
	ImmuClone Anti-D rapid	?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
	NOVACLONE	0	0	30

NEO Iris Operator Manual

Assay	Reaction	eaction Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	< =
	Anti-D	?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
ABD6CDE_I	ImmuClone Rh-Hr Control	0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		? 30	58	
	ImmuClone Anti-A	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-B	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
	ImmuClone Anti-AB	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
ImmuClone Anti-D rapid		0	0	30
		?	30	76
	1+ (not reported)	N/A	N/A	
	2+ (not reported)	N/A	N/A	
Assay	Reaction	Grade	Lower Limit	Upper Limit
-------	-----------------------	----------------------	-------------	-------------
		3+ (not reported)	N/A	<= N/A
		4+	76	100
	ImmuClone Anti-CDE	0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	< =
		0	0	30
		?	30	76
	NOVACLONE	1+ (not reported)	N/A	N/A
ABD6CDE_N	Diluent Control	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
	NOVACLONE	0	0	30
	Anti-A	?	30	58

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
		?	30	76
	NOVACLONE Anti-B	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	? 30	76
	NOVACLONE Anti-AB	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
	ImmuClone	0	0	30
	Anti-D rapid	?	30	76

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
	ImmuClone Anti-CDE	?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		Glade	>	<=
		0	0	30
		?	30	76
ABD6K_I ImmuClone Rh-Hr Control	ImmuClone	1+ (not reported)	N/A	N/A
	Rh-Hr Control	2+ (not reported) N/A	N/A	N/A
		3+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		4+	76	100
		0	0	30
		?	30	58
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-B	2+ (not reported)	N/A	N/A
		3+ (not reported)	(not prted) N/A N	N/A
		4+	76	100
		0	0	30
		?	30	76
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-AB	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	ImmuClone Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	23
		?	23	50
	Automated immuClone Anti-Kell	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	50	80
		4+	80	100

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	< =
	NOVACLONE Diluent	0	0	30
ABD6K_N Diluent Control		?	30	76
	1+ (not reported)	N/A	N/A	

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		Crade	>	<=
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	NOVACLONE Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
		?	30	76
	NOVACLONE Anti-B	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
	NOVACLONE	?	30	76
	Anti-AB	1+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
	ImmuClone Anti-D rapid	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	23
		?	23	50
	Automated immuClone Anti-Kell	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	50	80
		4+	80	100

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	< =
ABO_REV	A1-Cell	0	0	23

Assay	Peaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	A2-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
	B-Cell	?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
	O-Cell	?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100

Assay	Peaction	Grada	Lower Limit	Upper Limit	
Abbreviation	Reaction	Glade	>	< =	
		0	0	30	
		?	30	76	
	NOVACLONE	1+ (not reported)	N/A	N/A	
	Diluent Control	2+ (not reported)	N/A	N/A	
		3+ (not reported)	N/A	N/A	
		4+	76	100	
	NOVACLONE Anti-A	0	0	30	
		?	30	58	
ABODFULL			1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A	
		3+ (not reported)	N/A	N/A	
		4+	58	100	
		0	0	30	
			?	30	76
	NOVACLONE	1+ (not reported)	N/A	N/A	
	Anti-B	2+ (not reported)	N/A	N/A	
		3+ (not reported)	N/A	N/A	

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		1+	76	<= 100
				20
		0	0	30
		?	30	/6
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
	NOVACLONE Anti-D	?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	AT-Cell	2+	35	50
		3+	50	76
		4+	76	100
	B-Cell	0	0	23

NEO Iris Operator Manual

Assay	Peaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	U-Cell	2+	35	50
		3+	50	76
		4+	76	100

Assay	Peaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	< =
		0	0	30
		?	30	76
	ImmuClone Rh-Hr Control	1+ (not reported)	N/A	N/A
ABODFULL2		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
	ImmuClone	0	0	30
	Anti-A	?	30	58

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		1+ (not reported)	N/A	<= N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-B	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	ImmuClone Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
	NOVACLONE	0	0	30
	Anti-D	?	30	76

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
	A1-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
	P. Call	1+	28	35
	D-Cell	2+	35	50
		3+	50	76
		4+	76	100
		0	0	23
		?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	100

## Assay Reagent Component Grid

Reagents and Microplates	ABDLONG	ABDLONG2	ABDLONG3	ABDLONG4	ABDLONGK	ABDLONG2K	ABD6_I	ABD6_N	ABD6CDE_I	ABD6CDE_N	ABD6K_I	ABD6K_N	ABO_Rev	ABODFULL	ABODFULL2
Untreated Microplates (barcoded)	х	x	х	х	x	х	х	х	х	х	х	х	х	х	х
Reagent Red Blood cells (Referencells – Group A1)	х	x	x	х	x	х							х	х	х
Reagent Red Blood cells (Referencells – Group A2)	х	x	х	х	x	х							х		
Reagent Red Blood cells (Referencells – Group B)	х	x	х	х	x	х							х	х	х
Reagent Red Blood cells (Referencells – Group O)	х	x	x	х	x	х							х	х	х
NOVACLONE Anti-A	х		х		х			х		х		х		х	
NOVACLONE Anti-B	х		х		х			х		х		х		х	
NOVACLONE Anti-AB	х		х		х			х		х		х			
NOVACLONE Anti-D	Х	х	х	х	х	х	Х	Х						Х	Х
ImmuClone Anti-A		х		х		х	х		х		х				х
ImmuClone Anti-B		х		х		х	х		х		х				х

#### **Procedural Steps**

Reagents and Microplates	ABDLONG	ABDLONG2	ABDLONG3	ABDLONG4	ABDLONGK	<b>ABDLONG2K</b>	ABD6_I	ABD6_N	ABD6CDE_I	ABD6CDE_N	ABD6K_I	ABD6K_N	ABO_Rev	ABODFULL	ABODFULL2
ImmuClone Anti-AB		х		х		х	х		х		х				
ImmuClone Anti-D rapid	х	х	х	х	х	х	х	х	х	х	х	х		х	х
ImmuClone Rh-Hr Control	х	х	х	х	х	х	х		х		х				х
NOVACLONE Diluent Control								х		х		х		х	
ImmuClone Anti-CDE			х	х					х	х					
Automated immuClone Anti-Kell					х	х					х	х			
Autocontrol	х	х	х	х	х	х									

## Procedural Steps

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
ABDLONG	
	1. Bring all reagents and blood samples to 18–30°C before testing.
ABDLONG2	2. Centrifuge the blood samples to separate the plasma from the red blood
ABDI ONG3	cells and then remove the caps from those tubes. Process the donor
	segment blood samples, but do not centrifuge those samples.
ABDLONG4	3. Remove reagent vial caps.
	4. If applicable, add one stirball to each new vial of Referencells A1, B and O to
ABDLONGK	be used. Gently agitate each vial to resuspend the red blood cells.
ABDLONG2K	5. Load reagents, microplates, and blood samples, onto the NEO ins following
	6 Assign the appropriate assay to the blood samples either manually or
ABD6_I	following the upload worklist procedure
ABD6 N	7. Start the appropriate assay following the procedures in Chapter 6 –
1.000_11	Instrument Testing Operation. The NEO Iris automatically performs the
ABD6CDE_I	assays, and records and interprets blood sample results.
	8. At the completion of the NEO Iris assay, click the Results button on the main
ABD6CDE_N	menu bar to access the blood samples results.
ABD6K_I	
ABD6K_N	
ABO REV	
ABODFULL	
ABODFULL2	

## Results and Interpretations

### Introduction

The NEO Iris<sup>™</sup> generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported as Negative, Positive, Equivocal or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

## **Possible Test Well Results**

Possible test well results for assays include:

Indicator	Name	Description
+	Positive	The reaction value is greater than the cutoff value.
-	Negative	The reaction value is less than or equal to the cutoff value.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

## **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Ρ	ossible Test Interpretations
ABDLONG	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD / Mixed field ?, *INV*
ABDLONG2		Rh Interpretations:	RH+, RH-, NTD / Mixed field ?, *INV*
ABODFULL			
ABODFULL2			
ABDLONG3	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD / Mixed field ?, *INV*
ABDLONG4		Rh/CDE Interpretations	RH+ CDE+, RH- CDE-, RH- CDE+, NTD /
		Mixed field ?, *INV*	
ABDLONGK	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD/Mixed field ?, *INV*
ABDLONG2K		Rh Interpretations:	RH+, RH-, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*
ABD6_I	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD/Mixed field ?, *INV*
ABD6_N		Rh Interpretations:	RH+, RH-, NTD / Mixed field ?, *INV*
ABD6CDE_I	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD/Mixed field ?, *INV*
ABD6CDE_N		Rh/CDE Interpretations	RH+ CDE+, RH- CDE-, RH- CDE+, NTD / Mixed
		field ?, *INV*	
ABD6K_I	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD/Mixed field ?, *INV*
ABD6K_N		Rh Interpretations:	RH+, RH-, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*

Allaciment VI. NEO mis Operator Manua	Attachment	VI: N	EO Iris	Operator	Manual
---------------------------------------	------------	-------	---------	----------	--------

ABO_REV	+, -, ?, X	ABO Interpretations:	A, B, AB, NTD/Mixed field ?, *INV*
---------	------------	----------------------	------------------------------------

### In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results and interpretations used for the S3\_Cell assay.

ATTACHMENT VI: NEO IRIS OPERATOR MANUAL	VII-1
Copyrights and Disclaimers	VII-2
How this Attachment is Organized	VII-4
Assay Description	VII-6
Cutoffs and Reagent Components	VII-5
Procedural Steps	VII-6
Results and Interpretations	VII-9

## **Copyrights and Disclaimers**

©2018, Immucor. Inc. NEO Iris™ is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment VI for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrolde@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged.

No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

## In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

## **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (VII) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (AVII) in parentheses. The first three characters (NEO Iris) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

## Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

### Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
	Warning	Potentially damaging or dangerous outcomes if certain critical procedural steps are ignored or incorrectly executed.
īj	Consult instructions for use	

## Assay Description

### Introduction

### Intended Use

The S3\_Cell is an automated assay used for the detection of unexpected IgG antibodies to red blood cells for samples with low or no amount of albumin. This assay has to run in addition to the quality control assay QC3\_Cell before routine testing of blood specimens and can only function as an add-on test to the QC. <u>Sample addition is</u> <u>NOT verified when using the S3\_Cell assay</u>.

#### List of assays

The following table lists the described assay on the NEO Iris in a vertical assay configuration.

Assay Description	Assay short name	Used reagents	Microplates Used
	Capture-R Ready Screen		
RBC antibody screen	S3_Cell	<ol> <li>Capture LISS</li> <li>Capture-R Ready Indicator Red Cells</li> </ol>	Capture-R Ready Screen 3

### Limitations, Warnings and Notes

lcon	Description
	<u>Limitation:</u> Sample addition is NOT verified when using the S3_Cell assay as the Range Check of the Well filled verification step is turned off. Therefore each negative result requires confirmation by additional methodology. The S3_Cell assay is NOT intended as a stand-alone QC assay and does not replace the QC3_Cell Assay, the QC3_Cell assay must be run as the instrument quality control assay. This assay is not intended to be used for routine blood sample testing.

All further limitations, warnings and notes that relate to the described assay are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings and in Attachment I of the NEO Iris instrument.

## **Cutoffs and Reagent Components**

### List of Cutoffs

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
S3_Cell	Screening wells	0	0	20
		?	20	30
		1+	30	45
		2+	45	65

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		3+	65	90
		4+	90	100
	IgG coated positive control well	0	0	30
		?	-	-
		1+	30	45
		2+	45	65
	3+	65	90	
		4+	90	100

## **Procedural Steps**

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert.

Assay Button	Drief Our annie of Assess Desse dural Otana
Appreviation	Brief Synopsis of Assay Procedural Steps
S3 Cell	<ol> <li>Bring all reagents and blood samples to 18–30°C before testing.</li> </ol>
	2. Centrifuge the blood samples to separate the blood samples from the red blood
	cells/clot. Remove the caps from the blood sample tubes.
	<ol> <li>Remove the Capture-R<sup>®</sup> Ready-Screen<sup>®</sup> (3) microplate frame and the desired</li> </ol>
	number of Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (3) strips from the pouch.
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be
	Lood reasonable management of the subject of the real blood certains.
	<ol> <li>Load reagents, microplates, and blood samples onto the NEO ins following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	<ol> <li>Assign the S3_Cell assay to the blood samples, either manually or following the upload worklist procedure.</li> </ol>
	8. Start the S3_Cell assay following the procedures in Chapter 6 – Instrument
	Testing Operation. The NEO Iris automatically performs the S3_Cell assay, and
	records and interprets the blood sample results.
	9. At the completion of the NEO Iris S3_Cell assay, click the Results button on the
	main menu bar to access the blood samples results.

## Results and Interpretations

### Introduction

The NEO Iris<sup>™</sup> generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported as Negative, Positive, Equivocal or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

## **Possible Test Well Results**

Possible test well results for the S3\_Cell assay includes:

Indicator	Name	Description
+	Positive	The reaction value is greater than the positive cutoff value.
-	Negative	The reaction value is less than or equal to the negative control value.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.
X	Invalid	The invalid symbol indicates an error status for a given well. An Invalid test well results in generated if the instrument detects a processing error or process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.)

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations
S3_Cell	+, -, ?, X	Positive, Negative, No_Int, Ctrl Fail, *INV* <u>Note:</u> Ctrl Fail indicates that the per sample positive control well (4 <sup>th</sup> well in S3_Cell) has failed. This control well must react at least 3+ in order for the sample test result to be valid. If the control well is less than 3+, the the interpretation given will be Ctrl Fail.

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results and interpretations used for the Phenotyping assays.

ATTACHMENT VIII: NEO IRIS OPERATOR MANUAL	VIII-1
Copyrights and Disclaimers	VIII-2
How this Attachment is Organized	VIII-4
Assay Descriptions	VIII-6
Cutoffs and Reagent Components	VIII-10
Procedural Steps	VIII-32
Results and Interpretations	VIII-33

## Copyrights and Disclaimers

©2018, Immucor. Inc. NEO Iris™ is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment VIII for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrolde@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

## In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

## **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (VIII) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (AVIII) in parentheses. The first three characters (NEO Iris) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

## Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

## Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
$\mathbf{\Lambda}$	Warning	Potentially damaging or dangerous
		outcomes if certain critical
		procedural steps are ignored or
		incorrectly executed.
1 i	Consult instructions for	
للمل	use	

## Assay Description

## Introduction

### Intended Use

The described assays are designed to provide phenotyping status on RhCE, D and Kell. Assays are based on direct agglutination technology for the NEO Iris<sup>™</sup> automated analyzer.

### List of assays

The following table lists the described assays on the NEO Iris. The assays exist in a horizontal and/or vertical configuration.

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
Rh Blood Group	RHFORM_D	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-C(2)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-c(2)</li> </ol>	Untreated microplates
		<ol> <li>5. immuClone Anti-E(1)</li> <li>6. immuClone Anti-E(2)</li> <li>7. immuClone Anti-e(1)</li> <li>8. immuClone Anti-e(2)</li> <li>9. immuClone Anti-Kell</li> <li>10. automated ImmuClone Anti-Kell</li> <li>11. immuClone Anti-D rapid</li> <li>12. immuClone Rh-Hr Control</li> </ol>	

Assay Description	Assay short name	Used reagents	Microplates Used
	RHFORM_Cw	1. immuClone Anti-C(1)	Untreated microplates
		2. immuClone Anti-C(2)	
		3. immuClone Anti-c(1)	
		4. immuClone Anti-c(2)	
		5. immuClone Anti-E(1)	
Rh Blood		6. immuClone Anti-E(2)	
Group		7. immuClone Anti-e(1)	
		8. immuClone Anti-e(2)	
		9. immuClone Anti-Kell	
		10. automated ImmuClone Anti-Kell	
		11. immuClone Anti-Cw	
		12. immuClone Rh-Hr Control	
	PHENO16_2	1. immuClone Anti-C(2)	Untreated microplates
		2. immuClone Anti-c(2)	
Rh Blood Group		3. immuClone Anti-E(2)	
		4. immuClone Anti-e(2)	
		5. immuClone Anti-Kell	
		6. immuClone Rh-Hr Control	
Rh Blood Group	PHENO16_3	1. immuClone Anti-C(2)	
		2. immuClone Anti-c(2)	Untreated microplates
		3. immuClone Anti-E(2)	
		4. immuClone Anti-e(2)	
		5. automated ImmuClone Anti-Kell	
		6. immuClone Rh-Hr Control	

### **Assay Description**

Assay Description	Assay short name	Used reagents	Microplates Used	
	PHENO16_4	1. immuClone Anti-C(1)		
		2. immuClone Anti-c(1)		
Rh Blood		3. immuClone Anti-E(1)		
Group		4. immuClone Anti-e(1)	Untreated micropiates	
		5. automated ImmuClone Anti-Kell		
		6. immuClone Rh-Hr Control		
		1. immuClone Anti-C(1)		
		2. immuClone Anti-c(1)	Untreated microplates	
Rh Blood	CcEe	3. immuClone Anti-E(1)		
Group		4. immuClone Anti-e(1)		
		5. immuClone Rh-Hr Control		
Rh Blood Group	CcEe2	1. immuClone Anti-C(2)		
		2. immuClone Anti-c(2)		
		3. immuClone Anti-E(2)	Untreated microplates	
		4. immuClone Anti-e(2)		
		5. immuClone Rh-Hr Control		
Rh Blood Group	CDE	1. immuClone Anti-CDE		
		2. immuClone Rh-Hr Control	Untreated microplates	
	Blood grouping	Vertical/8 wells per strip		
Rh Blood Group	PHENO12_2	1. immuClone Anti-C(2)		
		2. immuClone Anti-c(2)		
		3. immuClone Anti-E(2)	Untropted microplates	
		4. immuClone Anti-e(2)	ontreated micropiates	
		5. immuClone Anti-Kell		
		6. immuClone Rh-Hr Control		

Assay Description	Assay short name	Used reagents	Microplates Used
Rh Blood Group	PHENO12_3	<ol> <li>immuClone Anti-C(2)</li> <li>immuClone Anti-c(2)</li> <li>immuClone Anti-E(2)</li> <li>immuClone Anti-e(2)</li> <li>automated ImmuClone Anti-Kell</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
Rh Blood Group	CcEe	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
Rh Blood Group	CcEe2	<ol> <li>immuClone Anti-C(2)</li> <li>immuClone Anti-c(2)</li> <li>immuClone Anti-E(2)</li> <li>immuClone Anti-e(2)</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
Rh Blood Group	CDE	1. immuClone Anti-CDE 2. immuClone Rh-Hr Control	Untreated microplates

## Limitations, Warnings and Notes

All limitations, warnings and notes that relate to the described assays are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings and in Attachment I of the NEO Iris instrument. No additional limitations that relate to the described assays were identified.
## Cutoffs and Reagent Components

### List of Cutoffs

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		Giude	>	<=
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-C(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone Anti-C(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
RHFORM_D		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-c(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
	ImmuClone	?	23.0	75.0
	Anti-c(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
	ImmuClone Anti-E(2)	0	0.0	23.0
		?	23.0	75.0
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-e(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
	ImmuClone	?	23.0	75.0
	Anti-e(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
	immuClone Anti-Kell Automated ImmuClone Anti-Kell	?	23.0	35.0
		1+(not reported)	N/A	N/A
		2+	35.0	50.0
		3+	50.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	50.0
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	50.0	80.0
		4+	80.0	99.9
		0	0.0	30.0
		?	30.0	76.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-D rapid	2+(not reported)	N/A	N/A
		3+	76.0	76.0
		4+	76.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone Rh-Hr Control	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		Ciddo	>	<=
		4+	80.0	99.9

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	$\begin{array}{c ccccccccccccccccccccccccccccccccccc$	N/A	
	Anti-C(1)		N/A	
		3+	75.0	80.0
	1	4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
RHFORM_Cw	ImmuClone	1+(not reported)	N/A	N/A
	Anti-C(2)	2+(not reported) N/A	N/A	
		3+	23.0         75.0           N/A         N/A           N/A         N/A           75.0         80.0           80.0         99.9           0.0         23.0           23.0         75.0           N/A         N/A           N/A         N/A           10.0         23.0           23.0         75.0           N/A         N/A           N/A         N/A	80.0
		4+	80.0	99.9
		0	0.0	230
		?	23.0	75.0
	ImmuClone Anti-c(1)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-c(2)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
	0     0       ?     23       1+(not reported)     N,	0.0	23.0	
		?	23.0	75.0
		1+(not reported)	N/A	N/A
	Anti-E(1)	reported) 2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(2)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
	-	0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-e(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-e(2)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	35.0
	immuClone Anti-Kell	1+(not reported)	N/A	N/A
		2+	35.0	50.0
	3+	50.0	80.0	
		4+	80.0	99.9
	Automated	0	0.0	23.0
	ImmuClone	?	23.0	50.0

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Giude	>	<=
	Anti-Kell	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	50.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
	ImmuClone Anti-Cw	?	23.0	75.0
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	immuClone Rh-Hr Control	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
PHENO12_2	ImmuClone	0	0.0	23.0

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
	Anti-C(2)	?	23.0	75.0
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone Anti-c(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone Anti-E(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
	ImmuClone	0	0.0	23.0
	Anti-e(2)	?	23.0	75.0

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Gidde	>	<=
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
	immuClone Anti-Kell	?	23.0	35.0
		1+(not reported)	N/A	N/A
		2+	35.0	50.0
		3+	50.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	immuClone Rh-Hr Control	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit
PHENO12_3	ImmuClone	0	0.0	23.0
Anti-C(2)	?	23.0	75.0	

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-c(2)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
	ImmuClone Anti-E(2)	?	23.0	75.0
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
	4+	80.0	99.9	
		0	0.0	23.0
	ImmuClone	?	23.0	75.0
	Anti-e(2)	1+(not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	neuction		>	<=
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	50.0
	Automated	1+(not reported)	N/A	N/A
	ImmuClone Anti-Kell	2+(not reported)	N/A	N/A
		3+	50.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	immuClone Rh-Hr Control	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	< =
		0	0.0	23.0
PHENO16_2	ImmuClone	?	23.0	75.0
	Anti-C(2)	1+(not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Giude	>	<=
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-c(2)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(2)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone Anti-e(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	35.0
	immuClone	1+(not reported)	N/A	N/A
	Anti-Kell	2+	35.0	50.0
		3+	50.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	immuClone	1+(not reported)	N/A	N/A
	Rh-Hr Control	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	01000	>	<=
		0	0.0	23.0
		? 23.0	75.0	
PHENO16_3	ImmuClone	1+(not reported)	d) N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	< =
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-c(2)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone Anti-E(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone Anti-e(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
	Automated	0	0.0	23.0

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit
	ImmuClone	?	23.0	50.0
	Anti-Kell	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	50.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
	immuClone	?	23.0	75.0
		1+(not reported)	N/A	N/A
Rh-Hr Control	2+(not reported)	N/A	N/A	
		3+	75.0	80.0
		4+	80.0	99.9

Assay	Reaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	< =
		0	0.0	23.0
		?	23.0	75.0
PHENO16_4	ImmuClone	1+(not reported)	N/A	N/A
	Anti-C(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
	ImmuClone	0	0.0	23.0

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Cidde	>	< =
	Anti-c(1)	?	23.0	75.0
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-e(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
	Automated	0	0.0	23.0
	ImmuClone	?	23.0	50.0

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
	Anti-Kell	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	50.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
	immuClone Rh-Hr Control	?	23.0	75.0
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		0	0.0	23.0
		? 23.0 1+(not reported) N/A	23.0	75.0
	ImmuClone		N/A	N/A
CcEe	Anti-C(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
	ImmuClone	0	0.0	23.0
	Anti-c(1)	?	23.0	75.0

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	neuction	Grade	>	<=
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(1)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone Anti-e(1)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
	immuClone	?	23.0	75.0
	Rh-Hr Control	1+(not reported)	N/A	N/A

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-C(2)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
CcEe2	ImmuClone Anti-c(2)	0	0.0	23.0
		?	23.0	75.0
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
	ImmuClone	0	0.0	23.0
	Anti-E(2)	?	23.0	75.0

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-e(2)	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9
		0	0.0	23.0
		?	23.0	75.0
	immuClone	1+(not reported)	N/A	N/A
	Rh-Hr Control	2+(not reported)	N/A	N/A
		3+	75.0	80.0
		4+	80.0	99.9

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
CDE	immuClone	0	0.0	30.0

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
	Anti-CDE	?	30.0	76.0
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9
			0.0	30.0
		?	30.0	76.0
	immuClana	1+(not reported)	N/A	N/A
	Rh-Hr Control	2+(not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9

## Assay Reagent Component Grid

Reagents and Microplates	RHFORM_D	RHFORM_Cw	PHENO12_2	PHENO12_3	PHENO16_2	PHENO16_3	PHENO16_4	CcEe	CcEe2	CDE
Untreated Microplates (barcoded)	х	х	х	х	х	х	х	х	х	х
ImmuClone Anti-C(1)	х	х					х	х		
ImmuClone Anti-C(2)	х	х	х	х	х	х			х	
ImmuClone Anti-c(1)	х	х					х	х		
ImmuClone Anti-c(2)	х	х	х	х	х	х			х	
ImmuClone Anti-E(1)	х	х					х	х		
ImmuClone Anti-E(2)	х	х	х	х	х	х			х	
ImmuClone Anti-e(1)	х	х					х	х		
ImmuClone Anti-e(2)	х	х	х	х	х	х			х	
ImmuClone Anti-Kell	х	х	х		х					
Automated ImmuClone Anti-Kell	х	х		х		х	х			
ImmuClone Anti-Cw		х								
ImmuClone Anti-CDE										х
ImmuClone Anti-D rapid	х									
ImmuClone Rh-Hr Control	х	х	х	х	х	х	х	х	х	х

## Procedural Steps

#### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert.

Brief Synopsis of Assay Procedural Steps
1. Bring all reagents and blood samples to 18–30°C before testing.
2. Centrifuge the blood samples to separate the plasma from the red blood
cells and then remove the caps from those tubes. Process the donor
segment blood samples, but do not centrifuge those samples.
3. Remove reagent vial caps.
4. Load reagents, microplates, and blood samples, onto the NEO Iris following
the procedures in Chapter 6 – Instrument Testing Operation.
5. Assign the appropriate assay to the blood samples either manually or
following the upload worklist procedure.
6. Start the appropriate assay following the procedures in Chapter 6 –
Instrument Testing Operation. The NEO Iris automatically performs the
assays, and records and interprets blood sample results.
7. At the completion of the NEO Iris assay, click the Results button on the main
menu bar to access the blood samples results.

## Results and Interpretations

### Introduction

The NEO Iris<sup>™</sup> generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported as Negative, Positive, Equivocal or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

### **Possible Test Well Results**

Possible test well results for assays include:

Indicator	Name	Description
+	Positive	The reaction value is greater than the cutoff value.
-	Negative	The reaction value is less than or equal to the cutoff value.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater that the negative cutoff value or equal o or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	F	Possible Test Interpretations
RHFORM_D	+, -, ?, X	CDE Interpretations:	ccD.ee, CCD.ee, CcD.ee, ccD.Ee, CCD.Ee, CcD.Ee,
			ccD.EE, CcD.EE, CCD.EE, ccddee, CCddee,
			Ccddee, ccddEe, CCddEe, CcddEe, ccddEE,
			CcddEE, CCddEE, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*
RHFORM_Cw	+, -, ?, X	CE and Cw Interpretat	ions: cc ee Cw+, CC ee Cw+, Cc ee Cw+, cc
			Ee Cw+, CC Ee Cw+, Cc Ee Cw+, cc EE
			Cw+, Cc EE Cw+, CC EE Cw+, cc ee Cw-
			CC ee Cw-, Cc ee Cw-, ccEe Cw-, CC Ee
			Cw-, Cc Ee Cw-, cc EE Cw-, CcEE Cw-,
			CC EE Cw-, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*
RHFORM_Cw	+, -, ?, X	CE Interpretations:	cc ee, CC ee, Cc ee, cc Ee Cw+, CC Ee,
*Alternative			Cc Ee, cc EE, Cc EE, CC EE,
Aurora File			
			NTD / Mixed field ?, *INV*
		Cw Interpretations:	Cw+, Cw-, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*

PHENO12_2	+, -, ?, X	CE Interpretations:	CC EE,
PHENO12_3			Cc ee, cc ee, NTD / Mixed field ?, *INV*
PHENO16_2		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*
PHENO16_3			
PHENO16_4			
CcEe	+, -, ?, X	CE Interpretations:	CC EE,
CcEe2			Cc ee, ccee, NTD / Mixed field ?, *INV*
CDE	+, -, ?, X	CDE Interpretations:	CDE+, CDE-, NTD / Mixed field ?, *INV*

\*Alternative result interpretation depending on installed aurora file version

### In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results and interpretations used for the S2\_Cell and SAb\_ID assays.

ATTACHMENT IX: NEO IRIS OPERATOR MANUAL	IX-1
Copyrights and Disclaimers	IX-2
How this Attachment is Organized	IX-4
Assay Descriptions	IX-6
Cutoffs and Reagent Components	IX-8
Procedural Steps	IX-9
Results and Interpretations	IX-11

## Copyrights and Disclaimers

©2020, Immucor. Inc. NEO Iris™ is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment IX for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (<u>doccontrol-de@immucor.com</u>).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (VIX) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (AVIX) in parentheses. The first three characters (NEO Iris) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

### Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

### Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
	Warning	Potentially damaging or dangerous
		outcomes if certain critical
		procedural steps are ignored or
		incorrectly executed.
1 îi	Consult instructions for	
	use	

## Assay Description

#### Introduction

#### Intended Use

The S2\_Cell assay is an automated screen assay used for the detection of unexpected IgG antibodies to red blood cells. The SAb\_ID assay is used as an antibody panel to identify unexpected IgG antibodies. Both, the S2\_Cell assay and the SAb\_ID assay contain an internal Negative control and an internal Positive control (containing a mixture of antibodies). Therefore an additional daily QC is not required for those as it is built into the test system.

#### List of assays

The following table lists the described assay on the NEO Iris in a vertical assay configuration.

Assay Description	Assay short name	Used reagents	Microplates Used
	Capture-R Ready Screen		
RBC antibody screen	S2_Cell	<ol> <li>Capture LISS</li> <li>Capture-R Ready Indicator Red Cells</li> <li>Capture-R Positive Control Serum (Weak)</li> <li>Capture-R Negative Control Serum</li> </ol>	Capture-R Ready- Screen (I and II)
RBC antibody panel	SAb_ID	<ol> <li>Capture LISS</li> <li>Capture-R Ready Indicator Red Cells</li> <li>Capture-R Positive Control Serum (Weak)</li> <li>Capture-R Negative Control Serum</li> </ol>	Capture-R Ready-ID

### Limitations, Warnings and Notes

lcon	Description
	Limitation: Sample addition is NOT verified when using the S2_Cell or SAb_ID assay
	as the Range Check of the Well filled verification step is turned off. Therefore each
	negative result requires confirmation by additional methodology. The assays
	contain internal Positive and Negative controls that need to pass in order to
	validate the run of the blood samples. This assay is not intended to be used for
	routine blood sample testing.

All further limitations, warnings and notes that relate to the described assay are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings and in Attachment I of the NEO Iris instrument.

## Cutoffs and Reagent Components

### List of Cutoffs

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
	All	0	0	20
		?	20	30
		1+	30	45
S2_Cell		2+	45	65
		3+	65	90
		4+	90	100
	Sample	0	0	20
		?	20	30
		1+	30	45
		2+	45	65
		3+	65	90
		4+	90	100
		0	0	30
	Positive control	?	-	-
		1+	30	45
SAb_ID		2+	45	65
		3+	65	90
		4+	90	100
	Negative control	0	0	30
		?	-	-
		1+	30	45
		2+	45	65
		3+	65	90
		4+	90	100

## Procedural Steps

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
S2 Cell	1. Bring all reagents and blood samples to 18–30°C before testing.
_	2. Centrifuge the blood samples to separate the blood samples from the red blood
	cells/clot. Remove the caps from the blood sample tubes.
	3. Remove the Capture-R <sup>®</sup> Ready-Screen <sup>®</sup> (I and II) microplate frame and the desired
	number of Capture-R ${ m e}$ Ready-Screen ${ m e}$ (I and II) strips from the pouch.
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture- $R^{\ensuremath{\$}}$ Ready Indicator Red Cells to be
	used. Gently agitate each vial to resuspend the red blood cells.
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the
	procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the S2_Cell assay to the blood samples, either manually or following the
	upload worklist procedure.
	8. Start the S2_Cell assay following the procedures in Chapter 6 – Instrument Testing
	Operation. The NEO Iris automatically performs the S2_Cell assay, and records and
	interprets the blood sample results.
	9. At the completion of the NEO Iris S2_Cell assay, click the Results button on the
	main menu bar to access the blood sample results.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
SAb ID	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the blood samples from the red blood
	cells/clot. Remove the caps from the blood sample tubes.
	3. Remove the Capture-R <sup>®</sup> Ready-ID <sup>®</sup> microplate frame and the desired number of
	Capture-R <sup>®</sup> Ready-ID <sup>®</sup> strips from the pouch.
	4. Remove reagent vial caps.
	5. Add one stirball to each new vial of Capture-R <sup>®</sup> Ready Indicator Red Cells to be
	used. Gently agitate each vial to resuspend the red blood cells.
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the
	procedures in Chapter 6 – Instrument Testing Operation.
	7. Assign the SAb_ID assay to the blood samples, either manually or following the
	upload worklist procedure.
	8. Start the SAb_ID assay following the procedures in Chapter 6 – Instrument Testing
	Operation. The NEO Iris automatically performs the SAb_ID assay, and records the
	blood sample results.
	9. At the completion of the NEO Iris SAb_ID assay, click the Results button on the
	mainmenu bar to access the blood sample results. If any wells are reported positive,
	you must cross-reference this data with the lot-specific Capture-R $^{\ensuremath{\mathbb{R}}}$ Ready-ID $^{\ensuremath{\mathbb{R}}}$ Master
	List to determine the antibody identification (if any exists).
## Results and Interpretations

### Introduction

The NEO Iris<sup>™</sup> generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported as Negative, Positive, Equivocal or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

## **Possible Test Well Results**

Possible test well results for the S2\_Cell and SAb\_ID assay include:

Indicator	Name	Description
+	Positive	The reaction value is greater than the positive cutoff value.
-	Negative	The reaction value is less than or equal to the negative control value.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.
x	Invalid	The invalid symbol indicates an error status for a given well. An Invalid test well results in generated if the instrument detects a processing error or process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.)

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations
S2_Cell	+, -, ?*, X	Positive, Negative, No_Int, *INV*
SAb_ID	+, -, ?*, X	Positive, Negative, No_Int, *INV*, Ctrl Fail
		Note: Ctrl Fail indicates that the per sample positive control well ( $15$ th
		well in SAb_ID) has failed. This control well must react at least 3+ in
		order for the sample test result to be valid. If the control well is less
		than 3+, then the interpretation given will be Ctrl Fail.

## **Attachment X: NEO Iris Operator Manual**

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results, interpretations used for the Antigen screening assays.

ATTACHMENT X: NEO IRIS OPERATOR MANUAL	X-1
Copyrights and Disclaimers	X-2
How this Attachment is Organized	X-4
Antigen Screening Assay Description	X-5
Antigen Screening Assay Cutoffs	X-7
Antigen Screening Assay Reagent Component Grid	X-7
Antigen Screening Assay Procedural Steps	X-8
Antigen Screening Assay Results and Interpretations	X-12

## Copyrights and Disclaimers

©2020, Immucor, Inc. NEO Iris<sup>®</sup> is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment X for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrolde@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

## In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

## **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (X) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (A-X) in parentheses. The first set of characters (NEO Iris\_EU) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

## Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

## Use of Icons

lcon	Type of Warning	Related to
	Warning	Potentially damaging or dangerous outcomes if certain critical
	, i a i i i g	procedural steps are ignored or incorrectly executed.
<u>[]i</u>	Consult instructions for use	

The following icons appear in this attachment to alert you of warnings or limitations of use.

## Antigen Screening Assay Description

## Introduction

#### Intended Use

Antigen screening assays are intended for the detection of rare antigens on red blood cells. Three categories of screening assays are available :

- Agglutination based antigen screening assay that incubates at room temperature.
- Agglutination based antigen screening assay that incubates at 37°C.
- Solid phase based antigen screening assay.

The Antigen screening assays make use of antisera selected by the customer, and the instrument system performs the assay processing steps using a standard set of parameters. Performance Characteristics have not been established. Because the type of antisera used in the assay is not prescribed by Immucor, the user must establish that the performance level is appropriate for their intended application.

#### List of assays

	Assay Button		
Assay Description	Abbreviation	Microplates Used	
Agglutination based antigen screening assay that	mAgScrPT	Untroated microplates	
incubates at room temperature.	IIIAgSCIKI	Untreated micropiates	
Agglutination based antigen screening assay that	m A a C ar 27	Untracted microplates	
incubates at 37°C.	magscrsz	Untreated micropiates	
Solid phase based antigen screening assay.	pAgScrAHG	Capture-R <sup>®</sup> Select	

## Limitations, Warnings and Notes

lcon	Description
	Limitation: The mAgScrRT, mAgScr37 and pAgScrAHG assays use antisera selected by end users, therefore Immucor has not established the performance characteristics of the customer selected antisera with the Iris system.
	It is the user's responsibility to identify and validate the antisera for use in conjunction with these assay test methods. As the user selects the antisera, Immucor has not performed any verification or validation to establish the specific performance characteristics with the antisera. Therefore, the users must establish the performance level and ensure it is appropriate for their specific application of the assay. Unless the end user has validated that the performance level is suitable for their use, additional measures should be considered to establish results are valid.
	Limitation:
	For the <b>mAgScrRT</b> , <b>mAgScr37</b> and <b>pAgScrAHG</b> assays: The barcode number used on the antisera vial will appear as the antisera ID in reports/results rather than the product name (e.g. Fy <sup>a</sup> ). You can choose to manually enter the product ID (e.g. Fy <sup>a</sup> ) and lot rather than scanning in the barcode or create a barcode with that name.
^	Limitation:
	For the mAgScrRT, mAgScr37 and pAgScrAHG assays: You will need to aliquot the correct amount of antisera into a test tube as directed in the Assay Procedural Steps. If insufficient volume of antisera is loaded, then a portion of the tests ordered will be invalidated due to insufficient liquid.
	Limitation:
	The maximum number of samples that can be antigen typed per microplate is 92 for the pAgScrAHG assay and 94 for the mAgScrRT and mAgScr37assays.
	Limitation:
	DAT positive red blood cells will give false positive results when used with the pAgScrAHG assay.
	Limitation:
	Capture-R <sup>®</sup> Ready Indicator Red Cells should be used no more than 24 hours after a stir ball has been added to the vial. Vials of reagents, other than Indicator Red Cells, that have remained continuously on the NEO Iris for 72 hours (3 days) should be removed and replaced with fresh vials. Vials of reagents, other than Indicator Red Cells, that are removed from the NEO Iris when not in use and refrigerated can be used up to their expiration dates.
	Users must determine the on-board usage time and stability for the antisera they apply for this assay.

## Antigen Screening Assay Cutoffs

### List of Cutoffs

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		0	0	30
		?	30	35
mAgScrRT	A 11	1+	35	55
mAgScr37	All	2+	55	70
		3+	70	80
		4+	80	99
		0	0	20
		?	20	40
	A 11	1+	40	50
pagscrang	All	2+	50	72
		3+	72	90
		4+	90	99

## Antigen Screening Assay Reagent Component Grid

Reagents & Microplates	mAgScrRT	mAgScr37	pAgScrAHG
Rh-Control	Х	Х	
Anti-D Rapid (ImmuClone)	Х	Х	
corQC™ EXTEND Standard	Х	Х	Х
Untreated microplates	Х	Х	
Capture <sup>®</sup> LISS			Х
Capture-R <sup>®</sup> Ready Indicator Red Cells			Х
Capture-R <sup>®</sup> Positive Control Serum			Х
Capture-R <sup>®</sup> Negative Control Serum			Х
Capture-R <sup>®</sup> Select			Х

## Antigen Screening Assay Procedural Steps

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
	Qualification of Screening Material
	1. Select the desired screening antiserum for the antigen of interest.
	2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.
	3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.
	4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.
	The target dilution to be used for the agglutination based antigen screening assays ( <b>mAgScrRT</b> and
	mAgScr37) is the highest dilution that gives at least a 2+ reaction in tube. However, a lower dilution
	may be used if desired.
	Note: If commercial antiserum is to be used, refer to the package insert for instructions on using the
	reagent for tube testing. If non-reagent material is to be used, refer to the following steps as a guide
	for tube testing.
	5. If not using a commercially prepared red blood cell suspension, prepare a 2-5% suspension of red
	blood cells.
	6. Add 1 drop of each antiserum dilution to separate test tubes
	7. Add 1 drop of the red blood cell suspension to each tube
	8. Mix the contents of each tube and centrifuge for 15-45s at 900-1000 x g.
mAgScrRT	9. Gently agitate each tube to suspend the red blood cell buttons. Examine for agglutination.
	A serve Des sectores (Change
	Assay Procedural Steps
	Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the
	caps from those tubes. Process the donor segment blood samples, but do not centrifuge those
	samples.
	a value of 45ul of antisera for each reaction plus 150 ul of antisera required as dead volume in the
	test tube (e.g. if 20 samples are being screeped for one antigen using the <b>mAgScrRT</b> assay 1.05 ml
	of antisera (45ul x 20 + 150ul) would be required )
	4 Remove reagent vial caps
	5 Load the diluted screening antiserum test tube* reagents, microplates, and blood samples (and/or
	donor samples), onto the NEO following the procedures in <b>Chapter 6 – Instrument Testing</b>
	Operation.
	Note: The barcode number printed on the original antiserum vial, if scanned or typed into the ID
	field, will appear as the antiserum ID in the reports/results rather than the product name (e.g. Fyª).

Abbreviation         Brief Synopsis of Assay Procedural Steps           You can choose to manually enter the product ID (e.g. Fyª), or create and scan a barcode with that name ID, rather than scanning in the barcode printed on the vial label.         6. Assign the mAgScrRT assay to the screening antiserum. This will bring up the Crossmatch dialog box. Enter the ID of each sample to be screened by using the hand held scanner or by manual entry. The mAgScrRT assay can also be ordered following the Worklist procedure.           7. Start the mAgScrRT assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO automatically performs the mAgScrRT assay, and records and interprets the screening results.           8. At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening matiserum for the antigen of interest.           2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.           3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.           4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.
<ul> <li>You can choose to manually enter the product ID (e.g. Fy<sup>a</sup>), or create and scan a barcode with that name ID, rather than scanning in the barcode printed on the vial label.</li> <li>Assign the mAgScrRT assay to the screening antiserum. This will bring up the Crossmatch dialog box. Enter the ID of each sample to be screened by using the hand held scanner or by manual entry. The mAgScrRT assay can also be ordered following the Worklist procedure.</li> <li>Start the mAgScrRT assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO automatically performs the mAgScrRT assay, and records and interprets the screening results.</li> <li>At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening results.</li> <li>Qualification of Screening Material</li> <li>Select the desired screening antiserum for the antigen of interest.</li> <li>Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> </ul>
<ul> <li>name ID, rather than scanning in the barcode printed on the vial label.</li> <li>6. Assign the mAgScrRT assay to the screening antiserum. This will bring up the Crossmatch dialog box. Enter the ID of each sample to be screened by using the hand held scanner or by manual entry. The mAgScrRT assay can also be ordered following the Worklist procedure.</li> <li>7. Start the mAgScrRT assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO automatically performs the mAgScrRT assay, and records and interprets the screening results.</li> <li>8. At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening results.</li> <li>Qualification of Screening Material</li> <li>1. Select the desired screening antiserum for the antigen of interest.</li> <li>2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ul>
<ul> <li>6. Assign the mAgScrRT assay to the screening antiserum. This will bring up the Crossmatch dialog box. Enter the ID of each sample to be screened by using the hand held scanner or by manual entry. The mAgScrRT assay can also be ordered following the Worklist procedure.</li> <li>7. Start the mAgScrRT assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO automatically performs the mAgScrRT assay, and records and interprets the screening results.</li> <li>8. At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening results.</li> <li>Qualification of Screening Material</li> <li>1. Select the desired screening antiserum for the antigen of interest.</li> <li>2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ul>
<ul> <li>box. Enter the ID of each sample to be screened by using the hand held scanner or by manual entry. The mAgScrRT assay can also be ordered following the Worklist procedure.</li> <li>7. Start the mAgScrRT assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO automatically performs the mAgScrRT assay, and records and interprets the screening results.</li> <li>8. At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening results.</li> <li>Qualification of Screening Material</li> <li>Select the desired screening antiserum for the antigen of interest.</li> <li>Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ul>
<ul> <li>The mAgScrRT assay can also be ordered following the Worklist procedure.</li> <li>7. Start the mAgScrRT assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO automatically performs the mAgScrRT assay, and records and interprets the screening results.</li> <li>8. At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening results.</li> <li>Qualification of Screening Material</li> <li>1. Select the desired screening antiserum for the antigen of interest.</li> <li>2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ul>
<ul> <li>7. Start the mAgScrRT assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO automatically performs the mAgScrRT assay, and records and interprets the screening results.</li> <li>8. At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening results.</li> <li>Qualification of Screening Material</li> <li>1. Select the desired screening antiserum for the antigen of interest.</li> <li>2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ul>
The NEO automatically performs the mAgScrRT assay, and records and interprets the screening results.         8. At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening results.         Qualification of Screening Material         1. Select the desired screening antiserum for the antigen of interest.         2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.         3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.         4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.
<ul> <li>results.</li> <li>8. At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening results.</li> <li><u>Qualification of Screening Material</u></li> <li>1. Select the desired screening antiserum for the antigen of interest.</li> <li>2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ul>
<ul> <li>8. At the completion of the NEO mAgScrRT assay, press the Results button on the Main Menu Bar to access the screening results.</li> <li><u>Qualification of Screening Material</u></li> <li>1. Select the desired screening antiserum for the antigen of interest.</li> <li>2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ul>
access the screening results.         Qualification of Screening Material         1. Select the desired screening antiserum for the antigen of interest.         2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.         3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.         4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.
<ul> <li><u>Qualification of Screening Material</u></li> <li>Select the desired screening antiserum for the antigen of interest.</li> <li>Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ul>
<ol> <li>Select the desired screening antiserum for the antigen of interest.</li> <li>Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ol>
<ol> <li>Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.</li> <li>Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ol>
<ol> <li>Choose a red blood cell sample with heterozygous expression of the antigen of interest.</li> <li>Test this red blood cell sample against each of the antiserum dilutions using the test tube method.</li> </ol>
4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.
The target dilution to be used for the agglutination based antigen screening assays ( <b>mAgScrRT</b> and
mAgScr37) is the highest dilution that gives at least a 2+ reaction in tube. However, a lower dilution
may be used if desired.
Note: If commercial antiserum is to be used, refer to the package insert for instructions on using the
reagent for tube testing. If non-reagent material is to be used, refer to the following steps as a guide
for tube testing.
5. If not using a commercially prepared red blood cell suspension, prepare a 2-5% suspension of red
blood cells.
6. Add 1 drop of each antiserum dilution to separate test tubes.
7. Add 1 drop of the red blood cell suspension to each tube.
8. Mix the contents of each tube and centrifuge for 15-45s at 900-1000 x g.
9. Gently agitate each tube to suspend the red blood cell buttons. Examine for agglutination.
Assay Procedural Steps
1. Bring all reagents and blood samples to 18–30°C before testing.
2. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the
caps from those tubes. Process the donor segment blood samples, but do not centrifuge those
samples.
3. Aliquot the correct amount of diluted screening antisera into a test tube. The <b>mAgScr37</b> assay uses
a value of 45µl of antisera for each reaction plus 150 µl of antisera required as dead volume in the
test tube. (e.g. if 20 samples are being screened for one antigen using the <b>mAgScr37</b> assay, 1.05 mi
of antisera. (45µl x 20 + 150µl) would be required.)
4. Remove reagent vial caps.
5. Load the diluted screening antiserum test tube <sup>*</sup> , reagents, microplates, and blood samples (and/or
Chapter 6 – Instrument Testing
Uperation.
field will appear as the apticerum ID in the reports (results rather than the preduct rates (a.e. f. a)
Now can choose to manually enter the product ID (a.g. [13]) as create and created with thet

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
	name ID, rather than scanning in the barcode printed on the vial label.
	6. Assign the <b>mAgScr37</b> assay to the screening antiserum. This will bring up the Crossmatch dialog
	box. Enter the ID of each sample to be screened by using the hand held scanner or by manual entry.
	The mAgScr37 assay can also be ordered following the Worklist procedure.
	7. Start the <b>mAgScr37</b> assay following the procedures in <b>Chapter 6 – Instrument Testing Operation</b> .
	The NEO automatically performs the <b>mAgScr37</b> assay, and records and interprets the screening
	results.
	8. At the completion of the NEO mAgScr37 assay, press the Results button on the Main Menu Bar to
	access the screening results.
	Qualification of Screening Material
	1. Select the desired screening antiserum for the antigen of interest.
	2. Prepare serial dilutions of the antiserum in test tubes with 6% BSA (Bovine Serum Albumin) in PBS.
	3. Choose a red blood cell sample with heterozygous expression of the antigen of interest.
	4. Test this red blood cell sample against each of the antiserum dilutions using the test tube method.
	5. The target dilution to be used for the antiglobulin based antigen screening assays ( <b>pAgScrAHG</b> ) is
	the highest dilution that gives at least a weak positive (w+) reaction in tube. However, a lower
	dilution may be used if desired.
	Note: If commercial antiserum is to be used, refer to the package insert for instructions on using the
	reagent for tube testing. If non-reagent material is to be used, refer to the following steps as a guide
	for tube testing.
	6. If not using a commercially prepared red blood cell suspension, prepare a 2-5% suspension of red
	blood cells.
	7. Add 1 drop of each antiserum dilution to separate test tubes.
	8. Add 1 drop of the red blood cell suspension to each tube.
	9. Mix the contents of each tube and incubate at 36-38°C for 20 minutes.
	10. Wash the red blood cells in each tube at least three times with large volumes of saline. Decant
pAgScrAHG	thoroughly after each wash.
	11. Add Anti-Human Globulin in the amount specified by the manufacturer's package insert. Mix the
	contents of each tube thoroughly.
	12. Centrifuge for 15-45 seconds at 900-1000 x g, or as described in the manufacturer's package insert.
	13. Gently agitate each tube to suspend the red blood cell buttons. Examine for agglutination.
	Assay Procedural Steps
	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the
	caps from those tubes. Process the donor segment blood samples, but do not centrifuge those
	samples.
	3. Aliquot the correct amount of diluted screening antisera into a test tube. The <b>pAgScrAHG</b> assay
	uses a value of 55 $\mu$ l of antisera for each reaction plus 150 $\mu$ l of antisera required as dead volume in
	the test tube. (e.g. if 20 samples are being screened for one antigen using the <b>pAgScrAHG</b> assay,
	1.25 ml of antisera. (55µl x 20 + 150µl) would be required.)
	4. Remove reagent vial caps.
	5. Load the diluted screening antiserum test tube*, reagents, microplates, and blood samples (and/or
	donor samples), onto the NEO following the procedures in Chapter 6 – Instrument Testing

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
	Operation.
	Note: The barcode number printed on the original antiserum vial, if scanned or typed into the ID
	field, will appear as the antiserum ID in the reports/results rather than the product name (e.g. Fy <sup>a</sup> ).
	You can choose to manually enter the product ID (e.g. Fy <sup>a</sup> ), or create and scan a barcode with that
	name ID, rather than scanning in the barcode printed on the vial label.
	6. Assign the <b>pAgScrAHG</b> assay to the screening antiserum. This will bring up the Crossmatch dialog
	box. Enter the ID of each sample to be screened by using the hand held scanner or by manual entry.
	The <b>pAgScrAHG</b> assay can also be ordered following the Worklist procedure.
	7. Start the <b>pAgScrAHG</b> assay following the procedures in <b>Chapter 6 – Instrument Testing Operation</b> .
	The NEO automatically performs the <b>pAgScrAHG</b> assay, and records and interprets the screening
	results.
	8. At the completion of the NEO <b>pAgScrAHG</b> assay, press the <b>Results</b> button on the Main Menu Bar to
	access the screening results.

## Antigen Screening Assay Results and Interpretations

### Introduction

The NEO Iris generates a result for each well read by the instrument and an interpretation of the results. The Antigen screening assays have predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative, Invalid or Equivocal. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

### **Possible Test Well Results**

Possible test well results for Antigen Screening Assays include:

Well Result	Name	Description				
+	Positive	The reaction value can be considered as positive.				
-	Negative	The reaction value can be considered as negative.				
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.				
Х	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).				

## **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations					
mAgScrRT	+,-,?,X	Sample Interpretations:	Negative, Positive, No_Int, *INV*				
mAgScr37		Control Interpretations:	Negative, Positive, No_Int, *INV**				
pAgScrAHG							

## Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test Interpretation	Description
INV	Invalid
No_Int	No Interpretation

## **Attachment XI: NEO Iris Operator Manual**

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results, interpretations and interface data used for the Non-ABO Titration assays.

Attachment XI: NEO Iris Operator Manual	XI-1
In this attachment	XI-1
Copyrights and Disclaimers	XI-2
How this Attachment is Organized	XI-3
Non-ABO Titration Assay Descriptions	XI-4
Non-ABO Titration Cutoffs and Reagent Components	XI-7
Non-ABO Titration Procedural Steps	XI-9
Non-ABO Titration Results and Interpretations	XI-10
Performance Characteristics	XI-12
Interface Specification Information	XI-13

## **Copyrights and Disclaimers**

©2020, Immucor, Inc. NEO Iris<sup>®</sup> (NEO Iris<sup>®</sup>, hereinafter "NEO Iris") is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XI for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany.

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged.

No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

## In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (XI) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a string hyphenated format combined with a four character attachment alpha-numeric identifier (A-XI) in parentheses. The first set of characters (NEO Iris\_EU) identifies the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

## Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

## Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use and to indicate where additional information is provided.

lcon	Type of Warning	Related to
	Warning	Potentially damaging or dangerous outcomes if certain critical procedural steps are ignored or incorrectly executed.
i	Consult instructions for use	
	Note	

## Non-ABO Titration Assay Description

## Introduction

The assays are for the determination of IgG alloantibody titers from plasma or serum on the NEO Iris automated analyzer.

The clinical interpretation and clinical cutoffs will be established by the user.

The assays are designed for the determination of the IgG antibody titers using Capture<sup>®</sup> solid phase red cell adherence technology against an antigen present on the chosen cells used to perform the tests. The instrument will display the reaction strength of each well.

#### **Description of assays**

The following table lists the Non-ABO Titration assays available on the NEO Iris. The assays employ Capture-R<sup>®</sup> Select technology using the red cell products Panoscreen<sup>®</sup>/Panoscreen<sup>®</sup> C<sup>w</sup> and Panoscreen<sup>®</sup> EXTEND as antigen expressing cells. The assays determine IgG alloantibody titers ranging from 2 to 256.

Assay Button Abbreviation	Red Cell reagent used for the test	Assay group
	Panoscreen I	
T_IgG_FT	Panoscreen C <sup>w</sup> I	
T laG P2	Panoscreen II	Titration assays using Panoscreen/Panoscreen C <sup>w</sup>
	Panoscreen C <sup>w</sup> II	cells
T laG P3	Panoscreen III	
'_igO_i o	Panoscreen C <sup>w</sup> III	
T_lgG_E1	Panoscreen EXTEND I	
T_lgG_E2	Panoscreen EXTEND II	
T_lgG_E3	Panoscreen EXTEND III	Titration assays using Panoscreen EXTEND cells
T_lgG_E4	Panoscreen EXTEND IV	
T_lgG_E5	Panoscreen EXTEND V	

## Limitations, Warnings and Notes

lcon	Description
	Samples tested positive on an antibody screen or antibody identification method and determined as <2 or negative in a titration assay should be considered to have a titer of 1.
	Limitation:
	Capture-R <sup>®</sup> Select plates: Strips removed from pouches should be used within 16 hours.
	Capture-R <sup>®</sup> Ready Indicator Red Cells should be used no more than 24 hours after a stir ball has been added to the vial. Vials of reagents other than Indicator red cells, that have remained continuously on the NEO Iris for 72 hours (3 days) should be removed and replaced with fresh vials. Vials of reagents, other than Indicator red cells, that are removed from the NEO Iris when not in use and refrigerated can be used up to their expiry dates.
•	Warning: Inspect all reagents for the presence of foam before placing on the instrument when performing any assay.
	Do not vigorously agitate reagents. Shaking will produce foam in the vial that can cause the Liquid Level Detection (LLD) function of the pipetting system to erroneously aspirate foam and/or air rather than reagent. This will produce incorrect results or an error.
	<u>Warning</u> : If you are using two or more instruments, then the specific reagent vials for each instrument must be dedicated for use on that single instrument to ensure correct reagent volume tracking. If the actual reagent volume (less than the software numeric volume) is not sufficient for the number of tests scheduled, the instrument will produce invalid results and samples will need to be rescheduled for testing.
	Warning: If you do not add the stirballs to the Capture-R <sup>®</sup> Indicator Red Cells suspensions, the Panoscreen/Panoscreen C <sup>w</sup> Cells or Panoscreen EXTEND Cells, the results may be invalid or incorrect. Do not touch the stirballs. You should add them directly to the cellular reagent vials using the dispenser provided. Contamination of cellular reagents can occur if the stirballs are touched.
	per vial.
	<u>Warning</u> : In rare occasions a possible result for a sample tested on a titration assay is "check for inconsistent grading". The instrument will not provide a titer result. Inconsistent grading implies that one or more well grading result within the dilution series of a sample are higher than the previous well grading result despite higher dilution (refer to table "List of Cutoffs"). A sample that shows the grading failure needs to be repeated.
	Warning: To prevent contamination of specimens that are to be tested multiple times over an extended period of time it is recommended to aliquot and freeze samples instead of loading one aliquot serveral times on and off the instrument.

lcon	Description
	Limitation: Samples that exhibit excessive hemolysis or lipemia, or are icteric, should not be tested on the instrument. Samples that exhibit a hemolysis concentration of more than 108 mg/dl in IgG assays must not be tested on the instrument, because they may generate erroneous results. Color check in IgG assays can detect excessive hemolysis and invalids the result. Samples that exhibit a triglyceride concentration of more than 273 mg/dl in IgG assays must not be tested on the instrument, because they may generate erroneous results. Icteric samples (conjugated bilirubin) are tested until a concentration of 52.4 mg/dl in IgG assays without showing erroneous results. Icteric samples (unconjugated bilirubin) are tested until a concentration of 26.3 mg/dl in IgG assays without showing erroneous results.
	<u>Warning:</u> The antigen expression on red cells used for the assays can be heterozygous or homozygous depending on the chosen cell. As the number of bound IgG antibodies per red cell is expected to be higher for cells which express the relevant antigen homozygously, the titer of an antibody determined on a homozygous cell might be higher than the titer determined on a heterozygous cell due to the nature of the cells.
	Warning: The antigens of red cell donors are displayed on the Panoscreen/ Panoscreen C <sup>w</sup> Master List accompanying each lot and are listed the Package insert for Panoscreen EXTEND. The correct cells expressing the antigen of interest need to be chosen for the test.

This list of notes, warnings and limitations is not intended as a substitute for the detailed package insert of the reagents used with the Non-ABO Titration assays.

## Non-ABO Titration Cutoffs and Reagent Components

### List of Cutoffs

Assay Abbreviation	Grade	Lower Limit >	Upper Limit <=
T_lgG_P1	0	0.0	20.0
T_IgG_P2 T_IgG_P3	1+	20.0	50.0
T_lgG_E1 T_lgG_E2	2+	50.0	72.0
T_lgG_E3	3+	72.0	90.0
T_lgG_E5	4+	90.0	99.9

All non-ABO Titration assays include an on-plate run control using Capture-R control sera. Cutoffs for the control wells are as displayed below:

Control cutoffs	Grade	Lower Limit >	Upper Limit <=		
	0	0.0	17.5		
	1+	17.5	50.0		
Capture-R Weak Positive and	2+	50.0	72.0		
Negative Control	3+	72.0	90.0		
	4+	90.0	99.9		

## Assay Reagent Component Grid

Reagents and Microplates	T_lgG_P1	T_lgG_P2	T_lgG_P3	T_lgG_E1	T_lgG_E2	T_lgG_E3	T_lgG_E4	T_lgG_E5
Capture-R <sup>®</sup> Weak Positive Control Serum	Х	Х	Х	Х	Х	Х	Х	Х
Capture-R <sup>®</sup> Negative Control Serum	Х	Х	Х	Х	Х	Х	Х	Х
Capture-R <sup>®</sup> Select	Х	Х	Х	Х	Х	Х	Х	Х
Panoscreen <sup>®</sup> I / Panoscreen <sup>®</sup> C <sup>w</sup> I	Х							
Panoscreen <sup>®</sup> II / Panoscreen <sup>®</sup> C <sup>w</sup> II		Х						
Panoscreen <sup>®</sup> III / Panoscreen <sup>®</sup> C <sup>w</sup> III			Х					
Panoscreen <sup>®</sup> Extend I				Х				
Panoscreen <sup>®</sup> Extend II					Х			
Panoscreen <sup>®</sup> Extend III						Х		
Panoscreen <sup>®</sup> Extend IV							Х	
Panoscreen <sup>®</sup> Extend V								Х
Capture-R <sup>®</sup> Ready Indicator Red Cells	Х	Х	Х	Х	Х	Х	Х	Х
Capture® LISS	х	х	Х	Х	Х	Х	Х	х

## Non-ABO Titration Procedural Steps

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
	<ol> <li>Bring all reagents and blood samples to 18–30°C before testing.</li> </ol>
	2. Centrifuge the blood samples to separate the plasma from the red blood cells.
	Remove the caps from the blood sample tubes. (50µl per sample is used for
T laG P1	Non-ABO Titration Assays)
	3. Remove reagent vial caps.
T_lgG_P2	4. Select the red cell expressing the antigen of interest from the provided lot
_	Masterlist (for Panoscreen/ Panoscreen C <sup>w</sup> ) or Package insert (Panoscreen
T_lgG_P3	EXTEND).
	5. Add one stindari to each riew vial of Parloscreen/Parlo
I_IgG_E1	Red Cells, Gently agitate each vial to resuspend the red blood cells
	6. Load reagents, microplates, and blood samples onto the NEO Iris following the
1_190_22	1. procedures in Chapter 6 – Instrument Testing Operation.
T IgG E3	7. Assign the appropriate assay(s) to the blood samples, either manually or
-0 -	following the upload worklist procedure.
T_lgG_E4	8. Start the appropriate assay(s) following the procedures in Chapter 6 -
	Instrument Testing Operation. The NEO Iris automatically performs the
T_lgG_E5	assay(s), and records and interprets blood sample results.
	9. Upon completion of the NEO Iris Non-ABO Titrations assay(s), click the Results
	button on the main menu bar to access the blood sample results.
	10. Assay endpoint: This is determined automatically as the last result above the
	cutoff, which is set at 20 for Non-ABO titration assays.

## Non-ABO Titration Results and Interpretations

## Introduction

The NEO Iris generates a result for each well read by the instrument and provides an interpretation of the results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

### **Possible Test Well Results**

Possible test well results for Non-ABO Titration assay includes:

Indicator	Name	Description
+	Positive	The reaction value is greater than the cutoff value.
-	Negative	The reaction value is less than or equal to the cutoff value.
Х	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

The Non-ABO Titration assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

### **Result Interpretations**

The result interpretation is based on the reaction pattern of test well results. Final titer interpretation is the same for each Non-ABO Titration assay  $(T_IgG_E1, ...)$  and is listed in the table below.

Non-ABO Titration assay		Dilution (per well)							
		2	4	8	16	32	64	128	≥256
	<2 or negative	-	-	-	-	-	-	-	-
	2	+	-	-	-	-	-	-	-
	4	+	+	-	-	-	-	-	-
	8	+	+	+	-	-	-	-	-
	16	+	+	+	+	-	-	-	-
Titer	32	+	+	+	+	+	-	-	-
	64	+	+	+	+	+	+	-	-
	128	+	+	+	+	+	+	+	-
	≥256	+	+	+	+	+	+	+	+

## **Performance Characteristics**

#### **Specific Performance Characteristics**

Method comparison studies were performed at two (2) external sites and at Immucor, Inc. as an internal site. Specimens were tested on NEO Iris. Test results were compared for agreement between the automated titer assays and results for manually prepared dilutions.

Comparison Non-ABO Titrations vs. Manual Doubling Dilutions	Equal or within ±1 Doubling Dilution			Equal or within ±2 Doubling Dilutions			
N	n	Agreement (%)	LCI* (%)	n	Agreement (%)	LCI* (%)	
66	62	93.9	86.7	66	100	95.6	

\* Agreement at the 95% one-sided lower confidence interval

#### Reproducibility

The reproducibility of the Non-ABO Titration Assays was evaluated at two (2) external sites and at Immucor, Inc. as an internal site. Each site tested three (3) samples per assay, representing low, medium and high titers covering the range of the assay. The samples were tested in triplicate per run, two (2) runs per day, for five (5) nonconsecutive days.

Reproducibility Summary of All Assay Results			Equal or wit ±1 Doubling D	thin ilution	Equal or within ±2 Doubling Dilutions			
Assay	N	n	Agreement (%)	LCI* (%)	n	Agreement (%)	LCI* (%)	
T_lgG_P1	90	90	100	95.9	90	100	95.9	
T_lgG_P2	90	90	100	95.9	90	100	95.9	
T_lgG_E2	90	90	100	95.9	90	100	95.9	
T_lgG_E5	90	90	100	95.9	90	100	95.9	

NEO Iris EU-001-100 (A-XI)

\* Agreement at the 95% one-sided lower confidence interval

## Interface Specification Information

Non-ABO Titration assays information for the interface specification, detailed below, is supplemental information to Attachment 2 for the NEO Iris Operator Manual. Refer to Attachment 2 for the NEO Iris Operator Manual for a full description of the interface specification.

Field	Description	Possible Values				
R 3.4	Assay Name	T_lgG_P1, T_lgG_P2, T_lgG_P3, T_lgG_E1, T_lgG_E2, T_lgG_E3, T_lgG_E4, T_lgG_E5				
R 4.1	Original graded result	-, 1, 2, 3, 4, X				
		Corresponding Assay	Possible Values			
R 4.2	Interpreted Result 1	T_lgG_P1 T_lgG_P2 T_lgG_P3 T_lgG_E1 T_lgG_E2 T_lgG_E3 T_lgG_E4 T_lgG_E5	<2 or negative, 2, 4, 8, 16, 32, 64, 128, ≥256, Check for inconsistent grading, *INV*			

### **Result Record and Field 4 Measurement Values**

### LIS Well Result Identification for Field 4 Measurement Values (R 4.1)

Assay	Test phase				
	1	Dilution 1:2			
T InG P2	2	Dilution 1:4			
T IgG P3	3	Dilution 1:8			
T_lgG_E1	4	Dilution 1:16			
T_lgG_E2	5	Dilution 1:32			
T_lgG_E3	6	Dilution 1:64			
T_lgG_E4	7	Dilution 1:128			
I_IgG_E5	8	Dilution 1:256			

# Attachment XIV: Capture-R Select based Screening Assays

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results and interpretations used for the Capture-R Select based Screening assays.

ATTACHMENT XIV: Attachment for Capture-R Select based	Screening
Assays	XIV-1
Copyrights and Disclaimers	XIV-2
How this Attachment is Organized	XIV-4
Description of Assays	XIV-6
Cutoffs and Reagent Components	XIV-12
Procedural Steps	XIV-14
Results and Interpretations	XIV-16
RedCell Labeling Sheet	XIV-18

## Copyrights and Disclaimers

©2021, Immucor. Inc. NEO Iris® is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XIV for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrolde@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged.

#### Attachment XIV: Capture-R Select based Screening Assays

No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

## In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

## **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (XIV) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (A-XIV) in parentheses. The first three characters (NEO Iris) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

## Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

## Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
	Warning	Potentially damaging or dangerous outcomes if certain critical procedural steps are ignored or incorrectly executed.
	Limitation	
	Note	
ĺĺ	Consult instructions for use	

## Description for Assays

### Introduction

#### Intended Use

The Capture-R Select based Screening Assays are intended for detection and identification of IgG red blood cell antibodies on the NEO Iris<sup>®</sup> and NEO<sup>®</sup> v2.0 automated systems.

#### **Principles of the Assays**

The described assays are designed to provide customers the option to select specific red cell reagents or characterized red cells from other sources, such as donor red blood cells, for automated 2 or 4 cell IgG antibody screening at 37°C to meet regional specific needs. The assays are intended to provide the customer the ability to select characterized red cells based on customer defined criteria. Since reagent red cells may come in different formulations, 2 and 4 cell antibody screening assays will be available for use of red cell suspensions with a concentration of 2-5% or with a lower concentration of 0.8% (+/-0.2%).

The IgG detection assays are based on Immucor's patented Capture-R<sup>®</sup> Select technology and are designed for use on the NEO Iris<sup>®</sup> and NEO<sup>®</sup> v2.0 automated analyzers.

## List of assays

TI f - II	states as the left of	1:-+-+	all a secolar a se		+ I		. TI			- · · · - · · · · · · · · · · · · · · ·	f: + :
	Mind tanie	lists the	descrined	accave	on the	N = () r = 0	: INA :	accavc e	avict in	a vertical	continuiration
	wing tubic	insts the	acscribea	ussuys .			). IIIC (	ussuys (		u vertieur	configuration
								,			5

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Vertical/8 wells per strip	
Antibody screening (2-5%)	SC_2C	<ol> <li>RedCell I*</li> <li>RedCell II*</li> <li>corQC<sup>®</sup> EXTEND Standard</li> <li>Capture<sup>®</sup> LISS</li> <li>Capture-R<sup>®</sup> Positive Control Serum (Weak)</li> <li>Capture-R<sup>®</sup> Negative Control Serum</li> <li>Capture-R<sup>®</sup> Ready Indicator Red Cells</li> </ol>	Capture-R <sup>®</sup> Select plates
Antibody screening (0.8%)	SC_2C_I	<ol> <li>RedCell 0.8 I*</li> <li>RedCell 0.8 II*</li> <li>corQC<sup>®</sup> EXTEND Standard</li> <li>Capture<sup>®</sup> LISS</li> <li>Capture-R<sup>®</sup> Positive Control Serum (Weak)</li> <li>Capture-R<sup>®</sup> Negative Control Serum</li> <li>Capture-R<sup>®</sup> Ready Indicator Red Cells</li> </ol>	Capture-R <sup>®</sup> Select plates

#### Attachment XIV: Capture-R Select based Screening Assays

Assay Description	Assay short name	Used reagents	Microplates Used
Antibody screening (2-5%)	SC_4C	<ol> <li>RedCell I*</li> <li>RedCell II*</li> <li>RedCell III*</li> <li>RedCell IV*</li> <li>corQC<sup>®</sup> EXTEND Standard</li> <li>Capture<sup>®</sup> LISS</li> <li>Capture-R<sup>®</sup> Positive Control Serum (Weak)</li> <li>Capture-R<sup>®</sup> Negative Control Serum</li> <li>Capture-R<sup>®</sup> Ready Indicator Red Cells</li> </ol>	Capture-R <sup>®</sup> Select plates
Antibody screening (0.8%)	SC_4C_I	<ol> <li>RedCell 0.8 I*</li> <li>RedCell 0.8 II*</li> <li>RedCell 0.8 III*</li> <li>RedCell 0.8 IV*</li> <li>corQC<sup>®</sup> EXTEND Standard</li> <li>Capture<sup>®</sup> LISS</li> <li>Capture-R<sup>®</sup> Positive Control Serum (Weak)</li> <li>Capture-R<sup>®</sup> Negative Control Serum</li> <li>Capture-R<sup>®</sup> Ready Indicator Red Cells</li> </ol>	Capture-R <sup>®</sup> Select plates

\*NOTE: The reagents described as "RedCell (0.8) I-IV" refer to the red cell reagents selected and validated by the customer for use in their laboratory. Immucor will provide 7-digit barcodes for each of the listed RedCell reagents. Those barcodes shall be attached to the selected red cells to enable their use on the NEO Iris<sup>®</sup> and NEO<sup>®</sup> v2.0 instruments. They will not contain any further information on lot numbers and expiry dates. Recording the relevant reagent information will be within the responsibility of the customer. A reagent tracking sheet provided with this Operator Manual Attachment shall facilitate simple tracking of used reagent lots and corresponding expiry dates.
## Limitations, Warnings and Notes

lcon	Description
	Warning: Red cell reagents suitable for the proposed test are required to be
	selected by the user. Select reagents based on the corresponding red cell
	Masterlist and further reagent information provided by the reagent's supplier and
	based on any applicable regulatory requirements. It is the users responsibility to
	select red cells with the phenotypes matching those specified in local regulations
	and guidelines applicable to red cell antibody screening techniques.
	Warning: Assays SC_2C and SC_4C are designed for the use of <i>RedCell</i> suspensions
	with a concentration of 2-5% while 0.8% (+/-0.2%) RedCell suspensions can be
	tested using the SC_2C_I and SC_4C_I assays. Deviation from the recommended
	concentrations may lead to false positive, false negative or invalid results.
	If red blood cells need to be diluted to meet the above requirements, dilute red
	blood cell concentrate in buffered saline.
	Red blood cells should be free of hemolysis. Fragmented RBC membranes will
	interfere with cell monolayer formation. If the red blood cell sample shows signs of
	degradation, i.e., slight to moderate (1+ - 2 +) hemolysis, wash sample red blood
	cells at least 2 times with buffered saline before preparing the desired red cell
	suspension. Samples exhibiting severe hemolysis ( $\geq$ 3+) cannot be used.
	Warning: On board time and reagent expiry of the RedCell reagents cannot be
	automatically tracked by the system as barcode labels do not contain the relevant
	information. The operator is required to track lot number, expiry date and first on
	board time/date for each <i>RedCell</i> vial used on the instrument. A <i>RedCell</i> labeling
	sheet is provided with this Operator Manual Attachment.
	Warning: Immucor does not provide an on board stability claim for any red cell
	reagents other than Immucor red cell products released for automated use on the
	NEO Iris <sup>®</sup> /NEO <sup>®</sup> v2.0. Immucor reagents for automated detection of IgG antibodies
	on NEO Iris <sup>®</sup> have a 72 hours on board stability claim. Use of Immucor reagent red
	cells is therefore recommended.

lcon	Description				
	Warning: Reagent fill volume for <i>RedCell</i> reagents is not automatically recognized				
	or tracked in the NEO software. When the <i>RedCell</i> reagent is loaded onto the				
	instrument, the operator is required to manually enter the correct fill volume in the				
	field next to the reagent name. Note that a dead volume of 1 mL needs to be				
	subtracted from the actual fill volume when entering the reagent information, i.e.				
	when the actual vial fill volume is 10 mL, manually type 9 mL in the corresponding				
	data field. To edit the fill volume, click on the relevant loading bay and reagent				
	rack and enter the correct fill volume minus the dead volume under Liquid				
	Properties Volume (see image below).				
	Identifiers and Liquid Properties				
	Back:         Prefix + Barcode Settings         Volume [ml]:         Information:				
	1 Reag1 0000421 0000 RedCell 0.81				
	2 Reag1 0000422				
	3 Reag1 0000423				
	4 Reagi 0000424				
•	Warring RodCallaward on the Sananing agains need to be visled in 10 mL slage				
	warning: <i>RedCells</i> used on the screening assays need to be vialed in 10 mL glass				
<u> </u>	Vials with the following specifications: 22 x 60 mm, thread 18 mm, flat bottom. Use				
	insufficient red cell volume or pipetter crashes				
•	Warning: Store all reagents used on the Canture-R Select based Screening Assays				
	according to the suppliers' specifications to prevent compromising of reagent				
<u> </u>	according to the suppliers' specifications to prevent compromising of reagent				
•	Warning: Do not nool red coll reagants from different vials or transfer red coll				
	warning. Do not pool red cell reagents non different vials of transfer red cell				
<b></b>	reagents from their original vial to prevent bacterial or chemical contamination.				
	Warning: In a case where both positive control wells obtain reaction scores below the positive cutoff the result for the Capture-R <sup>®</sup> Positive Control is "negative" or in				
	a case where both negative control wells obtain 3+ reactions the result for the				
	Capture-R® Negative Control is "positive" due to the definition in the DMS results				
	text. Although the failed control does not display the Ctrl_Fail notice, the run will				
	be invalidated by the system and no results interpretation will be provided.				
	Limitation: Capture-R Select based Screening Assays do only detect IgG antibodies				
	at a temperature of 37°C.				

#### Attachment XIV: Capture-R Select based Screening Assays

lcon	Description
	Note: It is recommended to qualify the use of the chosen <i>RedCells</i> by running
	a sample on the system with a known IgG antibody once per day.
	Note: The RedCell Labeling Sheet that is provided in this Operator Manual
	Attachment shall help the operator to track relevant RedCell reagent data and
	ensure correct labeling and usage of the red cell products

All further limitations, warnings and notes that relate to the described assay are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings and in Attachment I of the NEO Iris<sup>®</sup> instrument.

# Cutoffs and Reagent Components

## List of Cutoffs

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		0	0.0	17.5
		?	17.5	17.5
	Positive/ Negative	1+	17.5	50.0
	Control	2+	50.0	72.0
		3+	72.0	90.0
SC_2C		4+	90.0	99.9
SC_2C_1		0	0.0	20.0
		?	20.0	30.0
	RedCell (0.8) I/	1+	30.0	45.0
	RedCell (0.8) II	2+	45.0	65.0
		3+	65.0	90.0
		4+	90.0	99.9
		0	0.0	17.5
		?	17.5	17.5
	Positive/ Negative	1+	17.5	50.0
	Control	2+	50.0	72.0
		3+	72.0	90.0
SC_4C		4+	90.0	99.9
SC_4C_1		0	0.0	20.0
	RedCell (0.8) I/	?	20.0	30.0
	RedCell (0.8) II/	1+	30.0	45.0
	RedCell (0.8) III/	2+	45.0	65.0
	RedCell (0.8) IV	3+	65.0	90.0
		4+	90.0	99.9

### Assay Reagent Component Grid

Reagents and Microplates	sc_2C	sc_2c_I	SC_4C	SC_4C_I
Capture-R <sup>®</sup> Select plate	x	х	x	х
corQC <sup>®</sup> EXTEND Standard	x	x	x	x
Capture-R <sup>®</sup> Indicator Red Cells	x	x	x	x
Capture-R <sup>®</sup> Positive Control Serum (Weak)	x	x	х	x
Capture-R <sup>®</sup> Negative Control Serum	x	х	х	х
Capture <sup>®</sup> LISS	x	x	x	x
RedCell I*	x		x	
RedCell II*	x		x	
RedCell III*			х	
RedCell IV*			х	
RedCell 0.8 I*		х		х
RedCell 0.8 II*		x		х
RedCell 0.8 III*				x
RedCell 0.8 IV*				x

\*NOTE: The reagents described as "RedCell (0.8) I-IV" refer to the red cell reagents selected and validated by the customer for use in their laboratory.

# Procedural Steps

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements.

## Information on *RedCell* reagent labeling



Immucor provides barcode labels to enable the use of customer validated red cell reagents with concentrations of 2-5% or 0.8% for antibody screening purposes using the Capture-R Select based 2 and 4 Cell Screening assays. The following table lists the assignment of barcodes to the corresponding *RedCell* reagent.

RedCell reagent name	Barcode
RedCell I*	0000411
RedCell II*	0000412
RedCell III*	0000413
RedCell IV*	0000414
RedCell 0.8 I*	0000421
RedCell 0.8 II*	0000422
RedCell 0.8 III*	0000423
RedCell 0.8 IV*	0000424

To ensure that barcodes are correctly identified by the instrument, they shall only be used on reagent vials that align with the vial specifications described in the "Limitations, Warnings and Notes" section.

The listed barcodes do not automatically provide information on the lot, expiry dates and on board time of the *RedCell* reagent. For this purpose, a *RedCell* Labeling Sheet is attached to this Operator Manual Attachment. By filling in this sheet, tracking of lot numbers, expiry dates and further reagent information shall be ensured.

## Assay Procedural Steps

This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert.

Assay Button	Brief Synonsis of Assay Procedural Steps		
Abbreviation	Bhei Synopsis of Assay Procedural Steps		
	1. Bring all reagents and blood samples to 18–30°C before testing.		
	2. Select suitable red cell reagents to be used for antibody screening on the Capture-R		
	Select based Screening assays based on their antigen profile, red cell concentration		
	and suitability to run on the automated NEO Iris <sup>®</sup> / NEO <sup>®</sup> v2.0 platform.		
	If red cells are required to be diluted for use on the Screening assays, prepare 2-		
	5% or 0.8% (+/-0.2%) suspensions in buffered saline.		
	3. Attach barcode labels provided by Immucor to the chosen red cell reagent vials and		
	track lot number, expiry date and further reagent information using the RedCell		
	Labeling Sheet.		
	4. Centrifuge the blood samples to separate the plasma from the red blood cells.		
	Remove the caps from the blood sample tubes.		
	5. Remove the Capture-R $^{\ensuremath{\mathbb{R}}}$ Select microplate frame and the desired number of		
SC_2C	Capture-R <sup>®</sup> Select strips from the pouch.		
SC 2C I	6. Remove reagent vial caps.		
00_10_1	7. Add one stirball to each new vial of Capture- $R^{ extsf{8}}$ Ready Indicator Red Cells, corQC $^{ extsf{8}}$		
SC_4C	EXTEND cells and <i>RedCell</i> reagents to be used. Gently agitate each vial to		
	resuspend the red blood cells.		
SC_4C_I	8. Load reagents, microplates, and blood samples onto the NEO Iris following the		
	procedures in Chapter 6 – Instrument Testing Operation.		
	9. Manually enter the fill volumes of the <i>RedCell</i> reagents by clicking on the lane bay		
	on which the reagents are loaded and manually entering the fill volume minus the		
	vial deadvolume (1ml) under "Liquid Properties".		
	10. Assign the SC_2C, SC_2C_I, SC_4C and/ or SC_4C_I assay to the blood samples, either		
	manually or following the upload worklist procedure.		
	11. Start the SC_2C, SC_2C_I, SC_4C and/ or SC_4C_I assay following the procedures in		
	Chapter 6 – Instrument Testing Operation. The NEO Iris automatically performs the		
	SC_2C, SC_2C_I, SC_4C and/ or SC_4C_I assay, and records and interprets the blood		
	sample results.		
	12. At the completion of the NEO Iris SC_2C, SC_2C_I, SC_4C and/ or SC_4C_I assay, click		
	the Results button on the main menu bar to access the blood sample results.		

# Results and Interpretations

#### Introduction

The NEO Iris<sup>®</sup> generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris<sup>®</sup> assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported as Negative, Positive, Equivocal or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

#### **Possible Test Well Results**

Possible test well results for assays include:

Indicator	Name	Description
+	Positive	The reaction value is greater than the cutoff value.
-	Negative	The reaction value is less than or equal to the cutoff value.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

#### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations
SC_2C		
SC_2C_I	+ - 2 X	Run control interpretations: Positive, Negative, *INV*, Ctrl_Fail
SC_4C	τ, -, :, <b>Λ</b>	Result interpretations: Positive, Negative, No_Int, *INV*
SC_4C_I		

# RedCell Labeling Sheet

Original red cell reagent information			Information for use on NEO Iris			
Red cell reagent product name	Lot number	Cell concentration	Immucor barcode number (enter manually)	Instrument SN	Comments	Operator
	Expiry date	Vial fill volume	Fill volume entered in NEO software? (fill volume - 1ml)	First on board date and time		
			□ (<)			
Red cell reagent product name	Lot number	Cell concentration	Immucor barcode number (enter manually)	Instrument SN	Comments	Operator
	Expiry date	Vial fill volume	Fill volume entered in NEO software? (fill volume - 1ml)	First on board date and time		
			□ (✓)			
Red cell reagent product name	Lot number	Cell concentration	Immucor barcode number (enter manually)	Instrument SN	Comments	Operator
	Expiry date	Vial fill volume	Fill volume entered in NEO software? (fill volume - 1ml)	First on board date and time		
			□ (∽)			
Red cell reagent product name	Lot number	Cell concentration	Immucor barcode number (enter manually)	Instrument SN	Comments	Operator
	Expiry date	Vial fill volume	Fill volume entered in NEO software? (fill volume - 1ml)	First on board date and time		
			□ (∽)			

# **Attachment XV: NEO Iris Operator Manual**

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results and interpretations used for the ABO Rh assays.

Attachment XV: NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Description for ABODFULLH, ABORH_V and Baby_BG3 Assays	6
Cutoffs and Reagent Components	8
Procedural Steps	14
Results and Interpretations	16

# Copyrights and Disclaimers

©2021, Immucor. Inc. NEO Iris™ is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XIII for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrolde@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

# How this Attachment is Organized

#### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

#### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (XV) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (AXV) in parentheses. The first three characters (NEO Iris) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

### Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

## Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
$\mathbf{\Lambda}$	Warning	Potentially damaging or dangerous
		outcomes if certain critical
		procedural steps are ignored or
		incorrectly executed.
1 îi	Consult instructions for	
لتلمرا	use	

# Description for ABODFULLH, ABORH\_V and Baby\_BG3 Assays

#### Introduction

#### Intended Use

The described assays are designed to provide basic blood type and Rh status which can include both forward testing (of red blood cell) and reverse testing (of plasma/serum). Assays are based on direct agglutination technology for the NEO Iris<sup>™</sup> automated analyzer.

#### List of assays

The following table lists the described assays on the NEO Iris. The assays exist in a horizontal or vertical configuration.

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
		1. ImmuClone Rh-Hr Control 2. immuClone Anti-A 3. immuClone Anti-B	
Forward and Reverse ABORh Blood Grouping	ABODFULLH	<ol> <li>4. immuClone Anti-D rapid</li> <li>5. NOVACLONE Anti-D</li> <li>6. Referencells A1</li> </ol>	Untreated microplates
		7. Referencells B 8. Referencells O	

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Vertical/8 wells per strip	
Forward and Reverse ABORh Blood Grouping	ABORH_V	<ol> <li>ImmuClone Rh-Hr Control</li> <li>NOVACLONE Anti-A</li> <li>NOVACLONE Anti-B</li> <li>NOVACLONE Anti-A,B</li> <li>NOVACLONE Anti-D</li> <li>Referencells A1</li> <li>Referencells B</li> <li>Referencells O</li> </ol>	Untreated microplates
Forward ABORh Blood Grouping	Baby_BG3	<ol> <li>immuClone Rh-Hr Control</li> <li>NOVACLONE Anti-A</li> <li>NOVACLONE Anti-B</li> <li>NOVACLONE Anti-A,B</li> <li>NOVACLONE Anti-D</li> </ol>	Untreated microplates

## Limitations, Warnings and Notes

All limitations, warnings and notes that relate to the described assays are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings are in Attachment I of the NEO Iris instrument.

No additional limitations that relate to the described assays were identified.

# Cutoffs and Reagent Components

## List of Cutoffs

Assay	Reaction	Grade	Lower	Upper Limit
Abbreviation			Limit >	<=
		0	0.0	30.0
		? (not reported)	N/A	N/A
	ImmuClone Bh-Hr Control	1+	30.0	35.0
		2+	35.0	50.0
		3+	50.0	76.0
		4+	76.0	99.9
		0	0.0	30.0
		?	30.0	58.0
	immuClone Anti-A	1+ (not reported)	N/A	N/A
ABODFULLH		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58.0	99.9
	immuClone Anti-B	0	0.0	30.0
		?	30.0	76.0
		1+(not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9
		0	0.0	30.0
	immuClone	?	30.0	76.0
	Anti-D rapid	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9
		0	0.0	30.0
		?	30.0	76.0
	NOVACLONE Anti-D	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9
		0	0.0	23.0
	Referencells A1, B, O	?	23.0	28.0
		1+	28.0	35.0
		2+	35.0	50.0
		3+	50.0	76.0
		4+	76.0	99.9

Assay	Peaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	< =
		0	0.0	30.0
ABORH_V	ImmuClone	? (not reported)	N/A	N/A
	Rh-Hr Control	1+	30.0	35.0
		2+	35.0	50.0
		3+	50.0	76.0
		4+	76.0	99.9
	NOVACLONE	0	0.0	30.0

#### Attachment XV: NEO Iris Operator Manual

Assay	Prostion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
	Anti-A	?	30.0	58.0
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58.0	99.9
		0	0.0	30.0
		?	30.0	76.0
		1+ (not reported)	N/A	N/A
	Anti-B	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9
		0	0.0	30.0
		?	30.0	76.0
		1+ (not reported)	N/A	N/A
	Anti-A,B	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9
	NOVACLONE	0	0.0	30.0

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
	Anti-D	?	30.0	76.0
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9
		0	0.0	23.0
		?	23.0	28.0
Re	Referencells	1+	28.0	35.0
	A1, B, O	2+	35.0	50.0
		3+	50.0	76.0
		4+	76.0	99.9

Assay	Peaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	< =
		0	0.0	30.0
Baby_BG3	immuClone	? (not reported)	N/A	N/A
	Rh-Hr Control	1+	30.0	35.0
		2+	35.0	50.0
		3+	50.0	76.0
		4+	76.0	99.9
	NOVACLONE	0	0.0	30.0
	Anti-A	?	30.0	58.0

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58.0	99.9
		0	0.0	30.0
		?	30.0	76.0
	NOVACLONE Anti-B	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9
		0	0.0	30.0
	NOVACLONE Anti-A,B	?	30.0	76.0
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9
	NOVACLONE	0	0.0	30.0
	Anti-D	?	30.0	76.0

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	< =
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76.0	99.9

## Assay Reagent Component Grid

Reagents and Microplates	ABODFULLH	ABORH_V	Baby_BG3
Untreated Microplates (barcoded)	x	х	x
immuClone Rh-Hr Control	x	х	x
immuClone Anti-A	х		
immuClone Anti-B	х		
immuClone Anti-D rapid	x		
NOVACLONE Anti-A		х	x
NOVACLONE Anti-B		x	x
NOVACLONE Anti-A,B		х	x
NOVACLONE Anti-D	х	х	x
Referencells A1	x	x	
Referencells B	x	x	
Referencells O	x	х	

## Procedural Steps

## Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
ABODFULLH,	1. Bring all reagents and blood samples to 18–30°C before testing.
ABORH V	2. Centrifuge the blood samples to separate the plasma from the red blood
_	cells and then remove the caps from those tubes.
	3. Remove reagent vial caps.
	4. Add one stirball to each new vial of Referencells A1, B and O to be used.
	Gently agitate each vial to resuspend the red blood cells.
	5. Load reagents, microplates, and blood samples, onto the NEO Iris following
	the procedures in Chapter 6 – Instrument Testing Operation.
	6. Assign the appropriate assay to the blood samples either manually or
	following the upload worklist procedure.
	7. Start the appropriate assay following the procedures in Chapter 6 –
	Instrument Testing Operation. The NEO Iris automatically performs the
	assays, and records and interprets blood sample results.
	8. At the completion of the NEO Iris assay, click the Results button on the mair
	menu bar to access the blood samples results.

Assay Button			
Abbreviation	Brief Synopsis of Assay Procedural Steps		
Baby BG3	1.	Bring all reagents and blood samples to 18–30°C before testing.	
5-	2.	Centrifuge the blood samples to separate the plasma from the red blood	
	cells and then remove the caps from those tubes. Process the do		
		segment blood samples, but do not centrifuge those samples.	
	3.	Remove reagent vial caps.	
	4.	Load reagents, microplates, and blood samples, onto the NEO Iris following	
		the procedures in Chapter 6 – Instrument Testing Operation.	
	5.	Assign the appropriate assay to the blood samples either manually or	
		following the upload worklist procedure.	
	6.	Start the appropriate assay following the procedures in Chapter 6 –	
		Instrument Testing Operation. The NEO Iris automatically performs the	
		assays, and records and interprets blood sample results.	
	7.	At the completion of the NEO Iris assay, click the Results button on the main	
		menu bar to access the blood samples results.	

# Results and Interpretations

#### Introduction

The NEO Iris<sup>™</sup> generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported as Negative, Positive, Equivocal or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

#### **Possible Test Well Results**

Possible test well results for assays include:

Indicator	Name	Description		
+	Positive	The reaction value is greater than the cutoff value.		
-	Negative	The reaction value is less than or equal to the cutoff value.		
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater that the negative cutoff value or equal o or less than the positive cutoff value.		
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).		

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

#### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations		
ABODFULLH	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD / Mixed field ?, *INV*	
		Rh Interpretations:	RH+, RH-, NTD / Mixed field ?, *INV*	
ABORH_V	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD / Mixed field ?, *INV*	
		Rh Interpretations:	RH+, RH-, NTD / Mixed field ?, *INV*	
Baby_BG3	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD / Mixed field ?, *INV*	
		Rh Interpretations:	RH+, RH-, NTD / Mixed field ?, *INV*	

# **Attachment XVI: NEO Iris Operator Manual**

#### In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results, interpretations and interface data used for the Monoclonal assays for the detection of S, s, Fya, Fyb and k.

Attachment XVI: NEO Iris Operator Manual	. XVI-1
Copyrights and Disclaimers	. XVI-2
How this Attachment is Organized	. XVI-3
Description for monoclonal S/s, Fya, Fyb, k Assays	. XVI-4
Monoclonal S/s,Fy <sup>a</sup> ,Fy <sup>b</sup> ,k Assay Cutoffs and Reagent Components	sXVI-8
Monoclonal S/s, Fya, Fyb, k Assay Procedural Steps	XVI-12
Monoclonal S/s, Fya, Fyb, k Assay Results and Interpretations	XVI-13
Interface Specification Information	XVI-15
Performance Characteristics	XVI-16

# **Copyrights and Disclaimers**

©2022, Immucor, Inc. NEO Iris® is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XVI for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrolde@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged.

No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

#### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

#### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (XVI) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (A-XVI) in parentheses. The first set of characters (NEO Iris\_EU) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

#### Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

#### Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
	Warning	Potentially damaging or dangerous outcomes if certain critical procedural steps are ignored or incorrectly executed.
i	Consult instructions for use	

# Description for monoclonal S/s, Fya, Fyb, k Assays

#### Introduction

#### Intended Use

Anti-S (Monoclonal) Gamma-clone<sup>®</sup>, Anti-s (Monoclonal) Gamma-clone<sup>®</sup> and Anti-Fy<sup>b</sup> (Monoclonal) Gamma-clone<sup>®</sup> blood grouping reagents are intended for the detection of S, s, and Fy<sup>b</sup> antigens, respectively, on red blood cells by direct hemagglutination on the NEO Iris<sup>®</sup> automated analyzer.

Anti-Fy<sup>a</sup> (Monoclonal) Gamma-clone<sup>®</sup> and Anti-k (Monoclonal) Gamma-clone<sup>®</sup> blood grouping reagents are intended for the detection of Fy<sup>a</sup> and k antigens, respectively, on red blood cells by Capture-R<sup>®</sup> Select-based indirect agglutination on the NEO Iris<sup>®</sup> automated analyzer.

#### List of assays

The following table lists the monoclonal assays for the detection of S and s (AG\_Ssli\_m), Fy<sup>a</sup> (AG\_Fya\_m), Fy<sup>b</sup> (AG\_Fyb\_m) and k (AG\_k\_li\_m) on the NEO Iris. The AG\_Ssli\_m and AG\_Fyb\_m assays exist in a horizontal and in a vertical configuration.

|--|

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
Monoclonal S/s AG_Ssli_m		<ol> <li>Monoclonal Control</li> <li>Anti-S (Monoclonal) Gamma- clone<sup>®</sup></li> <li>Anti-s (Monoclonal) Gamma-clone<sup>®</sup></li> </ol>	Untreated microplates
Monoclonal Fy <sup>b</sup>	AG_Fyb_m	<ol> <li>Monoclonal Control</li> <li>Anti-Fy<sup>b</sup> (Monoclonal) Gamma- clone<sup>®</sup></li> </ol>	Untreated microplates
	Blood grouping	Vertical/8 wells per strip	
Monoclonal S/s	AG_Ssli_m	<ol> <li>Monoclonal Control</li> <li>Anti-S (Monoclonal) Gamma- clone<sup>®</sup></li> <li>Anti-s (Monoclonal) Gamma-clone<sup>®</sup></li> </ol>	Untreated microplates
Monoclonal Fy <sup>b</sup>	Monoclonal Fy <sup>b</sup> AG_Fyb_m 1. Monoclonal Control 2. Anti-Fy <sup>b</sup> (Monoclonal) Gamma- clone <sup>®</sup>		Untreated microplates
Monoclonal Fy <sup>a</sup>	AG_Fya_m	<ol> <li>Monoclonal Control</li> <li>Anti-Fy<sup>a</sup> (Monoclonal) Gamma- clone<sup>®</sup></li> </ol>	Capture-R <sup>®</sup> Select plates

Assay Description	Assay short name	Used reagents	Microplates Used
Monoclonal Cellano (k)	AG_k_li_m	<ol> <li>Monoclonal Control</li> <li>Anti-k (Monoclonal) Gamma-clone<sup>®</sup></li> </ol>	Capture-R <sup>®</sup> Select plates

## Limitations, Warnings and Notes

lcon	Description
	<u>Note</u> : The purpose of the Monoclonal Control for AG_Ssli_m, AG_Fya_m, AG_Fyb_m and AG_k_li_m assay is to serve as a sample control for Immucor low protein blood grouping reagents. It is expected to indicate those sample related conditions that could lead to spontaneous agglutination with low protein reagents and, therefore, a false positive interpretation for the test. When the Monoclonal Control well yields a positive result, then the instrument would not report the blood type.
	Limitation: Vials of reagents, that have remained continuously on the NEO Iris for 72 hours (3 days) should be removed and replaced with fresh vials. Vials of reagents that are removed from the NEO Iris when not in use and refrigerated can be used up to their expiration dates.
	Limitation: The NEO Iris cannot reliably detect hemagglutination reactions that are graded as 1+ or less in test tube methodology. The NEO Iris does not generate an interpretation of mixed-field. Such a mixed-field reaction will be interpreted as positive, negative, or equivocal.
	Limitation for AG_Fyb_m: The expression of certain red blood cell antigens may diminish in strength during storage, particularly in EDTA samples. Better results will be obtained with fresh samples.
	<u>Limitation</u> : The grading of reactions on the NEO Iris must only be regarded as an approximation when compared to off-line visual grading by laboratory technical staff.
	Limitation for AG_Ssli_m and AG_k_li_m: Samples that exhibit excessive hemolysis or lipemia, or are icteric, should not be tested on the instrument. Samples that exhibit a hemolysis concentration of more than 1034.4 mg/dl must not be tested on the assays, because they may generate erroneous results. Samples that exhibit a triglyceride concentration of more than 1536.0 mg/dl in the assays must not be tested on the instrument, because they may generate erroneous results. Icteric samples (conjugated bilirubin) are tested until a concentration of 41.6 mg/dl without showing erroneous results. Icteric samples (unconjugated bilirubin) are tested until a concentration of 40.0 mg/dl without showing erroneous results.

lcon	Description
	Limitation for AG Fya m and AG Fyb m: Samples that exhibit excessive hemolysis or lipemia, or are icteric, should not be tested on the instrument. Samples that exhibit a hemolysis concentration of more than 1008.0 mg/dl must not be tested on the assays, because they may generate erroneous results. Samples that exhibit a triglyceride concentration of more than 1794.0 mg/dl in the assays must not be tested on the instrument, because they may generate erroneous results. Icteric samples (conjugated bilirubin) are tested until a concentration of 44.8 mg/dl without showing erroneous results. Icteric samples (unconjugated bilirubin) are tested until a concentration of 40.8 mg/dl without showing erroneous results.
	Limitation for AG_Fyb_m: The Fyx phenotype may not be detected by the Immucor Anti-Fy <sup>b</sup> (Monoclonal) Gamma-clone reagent.

# Monoclonal S/s, Fy<sup>a</sup>, Fy<sup>b</sup>, k Assay Cutoffs and Reagent Components

#### List of Cutoffs

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		0	0	25
	Monoclonal control	? (not reported)	N/A	N/A
		1+	25	35
		2+	35	50
		3+	50	80
AG_Ssli_m		4+	80	99
	Anti-S	0	0	25
	(Monoclonal) Gamma-clone <sup>®</sup>	?	25	28
		1+	28	35
	Anti-s (Monoclonal) Gamma-clone <sup>®</sup>	2+	35	50
		3+	50	80
		4+	80	99

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		0	0	25
	Monoclonal control	? (not reported)	N/A	N/A
		1+	25	35
		2+	35	50
		3+	50	72
AG_Fyb_m		4+	72	99
	Anti-Fyb (Monoclonal)	0	0	25
		?	25	28
		1+	28	35
	Gamma-clone <sup>®</sup>	2+	35	50
		3+	50	72
		4+	72	99

Assay	Reaction	Grade	Lower Limit >	Upper Limit
AG_Fya_m	Positive and negative control	0	0	17.5
		? (not reported)	N/A	N/A
		1+	17.5	50
		2+	50	72
		3+	72	90
		4+	90	99
	Monoclonal Control	0	0	20
		? (not reported)	N/A	N/A
		1+	20	50
		2+	50	72
		3+	72	90
		4+	90	99
	Anti-Fya (Monoclonal) Gamma-clone®	0	0	20
		?	20	40
		1+	40	50
		2+	50	72
		3+	72	90
		4+	90	99

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
AG_k_li_m	Positive and negative control	0	0	17.5
		? (not reported)	N/A	N/A
		1+	17.5	50
		2+	50	72
		3+	72	90
		4+	90	99
	Monoclonal Control	0	0	20
		? (not reported)	N/A	N/A
		1+	20	50
		2+	50	72
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
-----------------------	--	-------	---------------	-------------------
		3+	72	90
		4+	90	99
	Anti-k (Monoclonal) Gamma-clone®	0	0	20
		?	20	40
		1+	40	50
		2+	50	72
		3+	72	90
		4+	90	99

## Assay Reagent Component Grid

Reagents and Microplates	AG_Ss_li_m	AG_Fyb_m	AG_Fya_m	AG_k_li_m
Untreated Microplates (barcoded)	х	х		
Capture-R <sup>®</sup> Select plates			x	х
Specimen Diluent	х	х		
Monoclonal Control	х	х	х	х
Anti-S IgM (Monoclonal) Gamma-clone®	х			
Anti-s IgM (Monoclonal) Gamma-clone®	х			
Anti-Fyb IgM (Monoclonal) Gamma-clone®		х		
Anti-Fya IgG (Monoclonal) Gamma-clone®			x	
Anti-k IgG (Monoclonal) Gamma-clone®				х
corQC <sup>®</sup> Extend Standard			x	х
Capture-R <sup>®</sup> Positive and Negative Control Sera			x	х
Capture <sup>®</sup> LISS			х	х
Capture-R <sup>®</sup> Ready Indicator Red Cells			x	х

## Monoclonal S/s, Fya, Fyb, k Assay Procedural Steps

### Before you begin

i

You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
AG_Ss_li_m, AG_Fyb_m	<ol> <li>Bring all reagents and blood samples to 18–30°C before testing.</li> <li>Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.</li> </ol>
	3. Remove reagent vial caps.
	<ol> <li>Load reagents, microplates, and blood samples, onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	<ol> <li>Assign the assay (AG_Ssli_m / AG_Fyb_m) to the blood samples either manually or following the upload worklist procedure.</li> </ol>
	<ol> <li>Start the AG_Ssli_m / AG_Fyb_m assay following the procedures in Chapter 6         <ul> <li>Instrument Testing Operation. The NEO Iris automatically performs the AG_Ssli_m / AG_Fyb_m assay, and records and interprets blood sample results.</li> </ul> </li> </ol>
	<ol> <li>At the completion of the NEO Iris assays, click the Results button on the main menu bar to access the blood samples results.</li> </ol>
AG_Fya_m, AG_k_li_m	<ol> <li>Bring all reagents and blood samples to 18–30°C before testing.</li> <li>Centrifuge the blood samples to separate the plasma from the red blood cells and then remove the caps from those tubes. Process the donor segment blood samples, but do not centrifuge those samples.</li> </ol>
	<ol> <li>Remove the Capture-R<sup>®</sup> Select microplate frame and the desired number of Capture-R<sup>®</sup> Select strips from the pouch.</li> </ol>
	<ol> <li>Remove reagent vial caps.</li> <li>Add one stirball to each new vial of Capture-R<sup>®</sup> Ready Indicator Red Cells or CorQC<sup>®</sup> Extend cells to be used. Gently agitate each vial to resuspend the red blood cells.</li> </ol>
	<ol> <li>Load reagents, microplates, and blood samples, onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> </ol>
	<ol> <li>Assign the assay (AG_Fya_m / AG_k_li_m) to the blood samples either manually or following the upload worklist procedure.</li> </ol>
	<ol> <li>Start the AG_Fya_m / AG_k_li_m assay following the procedures in Chapter 6         <ul> <li>Instrument Testing Operation. The NEO Iris automatically performs the AG_Fya_m / AG_k_li_m assay, and records and interprets blood sample results.</li> </ul> </li> </ol>
	<ol><li>At the completion of the NEO Iris assays, click the Results button on the main menu bar to access the blood samples results.</li></ol>

## Monoclonal S/s, Fya, Fyb, k Assay Results and Interpretations

### Introduction

The NEO Iris<sup>®</sup> generates a result for each well read by the instrument and an interpretation of the results. The Monoclonal S/s, Fya, Fyb, k Assays have predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative, Invalid or Equivocal. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

### Possible Test Well Results

Possible test well results for Monoclonal S/s, Fya, Fyb, k Assays include:

Indicator	Name	Description
+	Positive	The reaction value can be considered as positive.
-	Negative	The reaction value can be considered as negative.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

#### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations
AG_Ssli_m	Monoclonal Control: +(1, 2, 3, 4), -, X S/s: +(1, 2, 3, 4), -, ? <sup>1</sup> , X	<u>S/s interpretations:</u> S+, S-, s+, s-, NTD, Ctrl failure, *INV*
AG_Fyb_m	Monoclonal Control: +(1, 2, 3, 4), -, X Fyb: +(1, 2, 3, 4), -, ? <sup>1</sup> , X	<u>Fyb interpretations:</u> Fyb+, Fyb-, NTD, Ctrl failure, *INV*

Assay	Possible Well Results	Possible Test Interpretations
AG_Fya_m	Monoclonal Control: +(1, 2, 3, 4), -, X Fya: +(1, 2, 3, 4), -, ? <sup>1</sup> , X	<u>Fya interpretations:</u> Fya+, Fya-, NTD, Mono ctrl failure, *INV*
		Control interpretations: positive, negative, Ctrl failure, *INV*
AG_k_li_m	Monoclonal Control: +(1, 2, 3, 4), -, X k (Cellano): +(1, 2, 3, 4), -, ? <sup>1</sup> , X	<u>Cellano interpretations:</u> Cellano+, Cellano-, NTD, Mono ctrl failure, *INV*
		Control interpretations: positive, negative, Ctrl failure, *INV*

<sup>1</sup> Equivocal results of antibody wells may be edited. There cannot be an equivocal result for Monoclonal Control wells, wherefore negative control results are not editable.

### Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test Interpretation	Description
INV	Invalid
NTD	No Type Determined
Ctrl failure	Control failure
Mono ctrl failure	Negative control failure

## Interface Specification Information

The Monoclonal S/s, Fya, Fyb and k assay information for the interface specification, detailed below, is supplemental information to Attachment II for the NEO Iris Operator Manual (EU). Refer to Attachment II for the NEO Iris Operator Manual (EU) for a full description of the interface specification.

### **Result Record Field 4 Measurement Values**

Possible Values for Reaction F	Pattern		
-, 1, 2, 3, 4, ?, X			
Interpretation Values			
S Interpretation	s Interpretation		
S+, S-, NTD, Ctrl failure, *INV*	s+, s-, NTD, Ctrl failure, *INV*		
S+ s+, S- s+, S+ s-, S- s-, NTD, Ctrl failure, *INV*			
Fya+, Fya-, NTD, Mono ctrl failure, *INV*			
Fyb+, Fyb-, NTD, Ctrl failure, *INV*			
Cellano+, Cellano-, NTD, Mono	ctrl failure, *INV*		
	Possible Values for Reaction F-, 1, 2, 3, 4, ?, XInterpretation ValuesS InterpretationS+, S-, NTD, Ctrl failure, *INV*S+ s+, S- s+, S+ s-, S- s-, NTD,Fya+, Fya-, NTD, Mono ctrl failure, *INFyb+, Fyb-, NTD, Ctrl failure, *INCellano+, Cellano-, NTD, Mono ctrl		

<sup>1</sup> Applicable for alternative Aurora file with one Element

### LIS Well Result Identification for Result Field 4 Measurement Values

Assay	Те	est phase		
	1	Monoclonal Control		
AG_Ssli_m	2	Anti-S (Monoclonal) Gamma-clone®		
	3	Anti-s (Monoclonal) Gamma-clone®		
	1	Anti-Fyb (Monoclonal) Gamma-clone®		
AG_Fyb_III	2	Monoclonal Control		
	1	Anti-Fya (Monoclonal) Gamma-clone®		
AG_Fya_III	2	Monoclonal Control		
ACklim	1	Anti-k (Monoclonal) Gamma-clone <sup>®</sup>		
AG_K_II_III	2	Monoclonal Control		

## **Performance Characteristics**

### **Specific Performance Characteristics**

#### Method Comparison Studies:

Method comparison studies were performed at two (2) external sites and one (1) internal site. Immucor, Inc., as the manufacturer, was the internal site. The external sites were representative of blood collection establishments, hospital-based transfusion services, and/or clinical laboratories. Samples were tested with the reagent and also a comparator reagent. Test results were evaluated for agreement between reagents. Combined results from all sites are summarized in the following tables:

Note: Agreement between methods does not indicate which method is correct.

Assay performance represents NEO Iris and Galileo NEO (software version 3.1 or higher).

#### Anti-S (Monoclonal) Gamma-clone:

Initial Results N=975		Comparator Reagent			
		Positive	Negative		
	Decitive	<b>545</b>	2* 427	Positive Percent Agreement	99.81%
Anti-S	Positive	545		PPA (95% 1-Sided LCI)	99.13%
				Negative Percent Agreement	99.53%
	Negative	1**		NPA (95% 1-Sided LCI)	98.54%

Discordant samples were further genotyped by DNA molecular testing (PreciseType™ HEA BeadChip). \*One (1) sample gave mixed-field reaction with comparator reagent; one (1) sample gave mixed-field reaction on NEO Iris. \*\*Sample gave mixed-field reaction with comparator reagent. If these three (3) samples are excluded from the analysis, resolved PPA 99.58% (95% 1-sided LCI) and NPA 98.90% (95% 1-sided LCI).).

#### Anti-little s (Monoclonal) Gamma-clone:

Initial Results N=975		Comparat	tor Reagent		
		Positive	Negative		
	Depitive	700	0*	Positive Percent Agreement	100%
Anti-s	Positive	793	3"	PPA (95% 1-Sided LCI)	99.71%
			179	Negative Percent Agreement	98.35%
	Negative	0		NPA (95% 1-Sided LCI)	95.80%†

Discordant samples were further genotyped by DNA molecular testing (PreciseType™ HEA BeadChip). \*Two (2) samples initially typed s– with comparator reagent; repeat test with comparator reagent and DNA were s+. One (1) sample initially type s– with comparator reagent; repeat tests with comparator reagent and NEO Iris were s–, DNA was s+. Resolved NPA 99.44% (97.39%<sup>†</sup> 95% 1-sided LCI). <sup>†</sup>The Lower 99% CI was not met due to the lower number of s– samples in the population, N=179.

Anti-k	(Cellano)	(Monoclonal)	) Gamma-clone:
		-	

Initial Results N=834		Compara	tor Reagent		
		Positive	Negative		
	Desitive	010	0	Positive Percent Agreement	100%
Anti-k	Positive	819		PPA (95% 1-Sided LCI)	99.72%
				Negative Percent Agreement	100%
	Negative	Negative 0		NPA (95% 1-Sided LCI)	85.77%*

\*The Lower 99% CI was not met due to the lower number of k- samples in the population, N=15.

#### Anti-Fya (Monoclonal) Gamma-clone:

Initial Results N=974		Comparator Reagent			
		Positive	Negative		
				Positive Percent Agreement	100%
Anti-Fyª	Fositive	673	5	PPA (95% 1-Sided LCI)	99.66%
			Negative Percent Agreement	98.34%	
	Negative	0	296	NPA (95% 1-Sided LCI)	96.54%

Discordant samples were further genotyped by DNA molecular testing (PreciseType<sup>™</sup> HEA BeadChip). All five (5) samples initially typed Fy(a–) with comparator reagent; repeat test with comparator reagent and DNA were Fy(a+). Resolved NPA 100% (99.23% 95% 1-side LCI).

#### Anti-Fyb (Monoclonal) Gamma-clone:

Initial Results N=1236		Comparator Reagent			
		Positive	Negative		
	Desitive	604	7*	Positive Percent Agreement	98.89%
Anti-Fy <sup>b</sup>	FOSILIVE	624		PPA (95% 1-Sided LCI)	97.93%
				Negative Percent Agreement	98.84%
	Negative	Negative 7 <sup>†</sup>		NPA (95% 1-Sided LCI)	97.84%

Discordant samples were further genotyped by DNA molecular testing (PreciseType<sup>™</sup> HEA BeadChip). \*<sup>†</sup> one (1) sample in each category resolved in favor of the instrument result. Resolved PPA 99.05% (98.13% 95% 1-side LCI) and NPA 99.01% (98.05% 95% 1-side LCI).

<sup>†</sup>Three (3) samples typed as Fy(a+b+<sup>w</sup>) [Fy<sub>mod</sub> (Fy<sup>x</sup>)] due to the 265C>T SNP. Three (3) samples were unresolved as DNA testing gave low signal results (2) or was QNS for testing (1).

\*One (1) sample had debris in the well and tested Fy(b-) upon retest. Four (4) samples were falsely positive due to misaligned instrument ROIs (Region of Interest); three (3) samples were Fy(b-) upon retest, one (1) sample remained Fy(b+) upon retest. One (1) sample was unresolved as DNA testing gave low signal results.

#### **Precision Studies:**

#### Anti-S, Anti-s, Anti-k, Anti-Fya and Anti-Fyb (Monoclonal) Gamma-clone reagents:

Repeatability and Reproducibility were performed at three (3) sites by testing identical sample panels, containing positive and negative panel members in triplicate, testing two runs per day for five non-consecutive days. Results demonstrated 100% agreement for all positive and negative panel members.

# **Attachment XVIII: NEO Iris Operator Manual**

### In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results, interpretations and interface data used for the RHFORMEL and PHENO12 assays for the detection of phenotyping status on RhCE and Kell.

Attachment XVIII: NEO Iris Operator Manual	XVIII-1
Copyrights and Disclaimers	XVIII-2
How this Attachment is Organized	XVIII-3
Description for RHFORMEL and PHENO12	XVIII-4
Cutoffs and Reagent Components	XVIII-6
Procedural Steps X	VIII-11
Results and InterpretationsX	VIII-12
Interface Specification InformationX	VIII-14
Performance CharacteristicsX	VIII-15

## **Copyrights and Disclaimers**

©2022, Immucor, Inc. NEO Iris® is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XVIII for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrolde@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged.

No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (XVIII) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (A-XVIII) in parentheses. The first set of characters (NEO Iris\_EU) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

### Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

### Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

Icon	Type of Warning	Related to
	Warning	Potentially damaging or dangerous outcomes if certain critical procedural steps are ignored or incorrectly executed.
i	Consult instructions for use	

## **Description for RHFORMEL and PHENO12**

### Introduction

#### Intended Use

The described assays are designed to provide phenotyping status on RhCE and Kell. Assays are based on direct agglutination technology for the NEO Iris<sup>™</sup> automated analyzer.

#### List of assays

The following table lists the described assays on the NEO Iris. The assays exist in a horizontal or vertical configuration.

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
Rh Blood Group	RHFORMEL	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-C(2)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-c(2)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-E(2)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Anti-e(2)</li> <li>immuClone Anti-e(2)</li> <li>immuClone Anti-Kell (1)</li> <li>Automated ImmuClone Anti-Kell</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
	Blood grouping	Vertical/8 wells per strip	
Rh Blood Group	PHENO12	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Anti-Kell (1)</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates

### Limitations, Warnings and Notes

All limitations, warnings and notes that relate to the described assays are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings and in Attachment I of the NEO Iris instrument. No additional limitations that relate to the described assays were identified.

## **Cutoffs and Reagent Components**

### List of Cutoffs

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-C(1)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-C(2)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
	ImmuClone Anti-c(1)	0	0	23
		?	23	75
RHFORMEL		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-c(2)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-E(1)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-E(2)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-e(1)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
	ImmuClone Anti-e(2)	0	0	23
		?	23	75
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	35
	immuClone	1+ (not reported)	N/A	N/A
	Anti-Kell (1)	2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
		?	23	50
	Automated	1+ (not reported)	N/A	N/A
	Anti-Kell	2+ (not reported)	N/A	N/A
		3+	50	80
		4+	80	99
		0	0	23
	Rh-Hr Control	? (not reported)	N/A	N/A
		1+	23	35

#### Attachment XVIII: NEO Iris Operator Manual

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		2+	35	50
		3+	50	76
		4+	76	99

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-C(1)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone Anti-c(1)	1+ (not reported)	N/A	N/A
PHENO12		2+ (not reported)	N/A	N/A
-		3+	75	80
		4+	80	99
		0	0	23
	ImmuClone	?	23	75
		1+ (not reported)	N/A	N/A
	Anti-E(1)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
	ImmuClone	?	23	75
	Anti-e(1)	1+ (not reported)	N/A	N/A

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	35
	immuClone Anti-Kell (1)	1+(not reported)	N/A	N/A
		2+	35	50
		3+	50	80
		4+	80	99
	immuClone Rh-Hr Control	0	0	23
		? (not reported)	N/A	N/A
		1+	23	75
		2+	75	50
		3+	50	76
		4+	76	99

## Assay Reagent Component Grid

Reagents and Microplates	RHFORMEL	PHENO12
Untreated Microplates (barcoded)	х	х
ImmuClone Anti-C(1)	х	х
ImmuClone Anti-C(2)	х	
ImmuClone Anti-c(1)	х	х
ImmuClone Anti-c(2)	х	
ImmuClone Anti-E(1)	х	х
ImmuClone Anti-E(2)	х	
ImmuClone Anti-e(1)	х	х
ImmuClone Anti-e(2)	х	
ImmuClone Anti-Kell (1)	х	х
Automated ImmuClone Anti- Kell	х	
ImmuClone Rh-Hr Control	х	х

## **Procedural Steps**

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button	Brief Synonsis of Assay Procedural Stens
ADDIEVIATION	1 Drive all means and black and provide to 2000 before testing
RHFORMEL	1. Bring all reagents and blood samples to 18–30°C before testing.
	<ol><li>Centrifuge the blood samples to separate the plasma from the red blood cells</li></ol>
	and then remove the caps from those tubes. Process the donor segment blood
FILENOIZ	samples, but do not centrifuge those samples.
	3. Remove reagent vial caps.
	4. Load reagents, microplates, and blood samples, onto the NEO Iris following the
	procedures in Chapter 6 – Instrument Testing Operation.
	5. Assign the appropriate assay to the blood samples either manually or following
	the upload worklist procedure.
	6. Start the appropriate assay following the procedures in Chapter 6 – Instrument
	Testing Operation. The NEO Iris automatically performs the assays, and
	records and interprets blood sample results
	7 At the completion of the NEO Line account dick the Posults button on the main
	7. At the completion of the NEO ins assay, click the Results button on the main
	menu bar to access the blood samples results.

## Results and Interpretations

### Introduction

The NEO Iris<sup>®</sup>generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative, Invalid or Equivocal. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

### **Possible Test Well Results**

Possible test well results for assays include:

Indicator	Name	Description
+	Positive	The reaction value can be considered as positive.
-	Negative	The reaction value can be considered as negative.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations

Assay	Possible Well Results	Possible Test Interpretations
RHFORMEL PHENO12	ImmuClone Rh-Hr Control: +(1, 2, 3, 4), -, X Other reagents: +(1, 2, 3, 4), -, ? <sup>1</sup> , X	<u>CE Interpretations:</u> cc ee, CC ee, Cc ee, cc Ee, CC Ee, Cc Ee, cc EE, Cc EE, CC EE, NTD / Mixed field ?, NTD, *INV* <u>Kell Interpretations:</u> K+, K-, NTD / Mixed field ?, NTD, *INV*

<sup>1</sup> Equivocal results of antibody wells may be edited. There cannot be an equivocal result for ImmuClone Rh-Hr Control wells, wherefore negative control results are not editable.

### Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test Interpretation	Description
INV	Invalid
NTD	No Type Determined
NTD / Mixed field ?	No Type Determined / Mixed field ?
NTD NTD	No Type Determined No Type Determined

## Interface Specification Information

The RHFORMEL and PHENO12 assay information for the interface specification, detailed below, is supplemental information to Attachment II for the NEO Iris Operator Manual (EU). Refer to Attachment II for the NEO Iris Operator Manual (EU) for a full description of the interface specification.

### **Result Record Field 4 Measurement Values**

SUB-COMPONENT 1					
Assay	Possible Values for Reaction Pattern				
RHFORMEL PHENO12	-, 1, 2, 3, 4, ?, X				
SUB-COMPONENT 2					
Assay	Interpretation Values				
	CE Interpretation	Kell Interpretation			
RHFORMEL PHENO12	CC EE, Cc EE, cc EE, CC Ee, Cc Ee, cc Ee, CC ee, Cc ee, cc ee NTD, *INV*, NTD / Mixed field ?	K+, K-, NTD, *INV*, NTD / Mixed field ?			

### LIS Well Result Identification for Result Field 4 Measurement Values

Assay	Tes	Test phase				
	1	ImmuClone Anti-C(1)				
	2	ImmuClone Anti-C(2)				
	3	ImmuClone Anti-c(1)				
	4	ImmuClone Anti-c(2)				
	5	ImmuClone Anti-E(1)				
RHFORMEL	6	ImmuClone Anti-E(2)				
	7	ImmuClone Anti-e(1)				
	8	ImmuClone Anti-e(2)				
	9	ImmuClone Anti-Kell (1)				
	10	Automated Anti-Kell				
	11	ImmuClone Rh-Hr Control				
	1	ImmuClone Anti-C(1)				
	2	ImmuClone Anti-c(1)				
PHENO12	3	ImmuClone Anti-E(1)				
	4	ImmuClone Anti-e(1)				
	5	ImmuClone Anti-Kell (1)				
	6	ImmuClone Rh-Hr Control				

### **Performance Characteristics** Specific Performance Characteristics

#### Method Comparison Studies:

A method comparison study was performed at one (1) internal site at Immucor Medizinische Diagnostik GmbH. Result concordance of the updated RHFORMEL and PHENO12 assays (i.e. RHF\_PEO and PH12\_PEO) on the NEO Iris (methods under test) were compared to the validated CcEe assay (reference method 1) and AG\_K assay (reference method 2) on the NEO Iris.

To verify the result consistency of the method under test and the reference method the concordance of all test results was calculated by comparing the number of results in agreement with the total number of results. Minitab 20 was used for concordance calculations employing the 1-proportion exact method at the 95% lower bound confidence interval (LCI). The following tables display the result concordance.

Assay performance represents NEO Iris and Galileo NEO (software version 3.1 or higher).

#### Data analysis for method under test RHF\_PEO:

2x2 Contingency Table		Reference method		Analysis	
		Positive	Negative		
Method under test	Positive	34	0	Number of tests 258 Concordance 95% Lower Confidence interval 99.1%	* The PPA (95% 1-sided LCI) was 93.5% due to the low frequency of Kell positive antigens in the
	Negative	0	224	Negative percent agreement 100% (LCI: 93.5%*) Negative percent agreement 100% (LCI: 99.0%)	Caucausian population (9%). The percent agreement (point estimate) was 100%.

#### Data analysis for method under test PH12 PEO:

2x2 Contingency Table	Reference method		Analysis
	Positive	Negative	Number of tests

				253
	Positive	30	0	Concordance 95% Lower Confidence interval 99.1%
				Positive percent agreement
Method under test				
				Negative percent
				100% (LCI: 99.0%)
	Negative	0	223	

\* The PPA (95% 1-sided LCI) was 92.6% due to the low frequency of Kell positive antigens in the Caucausian population (9%). The percent agreement (point estimate) was 100%.

# Attachment XIX: Capture-R Select based Identification Assays

### In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results and interpretations used for the Capture-R Select based Identification assays.

A	TTACHMENT XIX: Attachment for Capture-R Select base	d Identification
A	ssays	XIX-1
	Copyrights and Disclaimers	XIX-2
	How this Attachment is Organized	XIX-4
	Description for 8-Cell ID Assays	XIX-6
	Cutoffs and Reagent Components	XIX-11
	Procedural Steps	XIX-14
	Results and Interpretations	XIX-17
	RedCell Labeling	
Sheet.	XIX-19	

## **Copyrights and Disclaimers**

©2022, Immucor. Inc. NEO Iris<sup>®</sup> is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XIX for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrolde@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged.

No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

## How this Attachment is Organized

### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (XIX) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (A-XIX) in parentheses. The first three characters (NEO Iris) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

### Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

### Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
	Warning	Potentially damaging or dangerous outcomes if certain critical procedural steps are ignored or incorrectly executed.
	Limitation	

	Note	
[]i	Consult instructions for use	

## Description for 8-Cell ID Assays

### Introduction

#### Intended Use

The Capture-R Select based Identification Assays are intended for the identification of unexpected IgG red blood cell antibodies on the NEO Iris automated systems.

#### **Principles of the Assays**

The assays are designed to provide customers the option to select specific red cell reagents for automated 8 cell IgG antibody identification at 37 °C to meet regional specific needs. The assays are intended to provide the customer the ability to select characterized red cells based on customer defined criteria. Since reagent red cells may come in different formulations, the 8 cell antibody identification assays will be available for use of red cell suspensions with a concentration of 2-5 % or with a lower concentration of 0.8 % (+/-0.2 %).

The IgG identification assays are based on Immucor's patented Capture-R<sup>®</sup> Select technology and are designed for use on the NEO Iris<sup>®</sup> and NEO<sup>®</sup> v2.0 automated analyzers.

#### List of assays

The following table lists the described assays on the NEO Iris. The assays exist in a vertical configuration.

Assay Description	Assay short name	Used reagents	Microplates Used
	Capture-R Select		
Antibody identification (2-5 %)	SC_8C	<ol> <li>RedCell III*</li> <li>RedCell IV*</li> <li>RedCell V*</li> <li>RedCell VI*</li> <li>RedCell VII*</li> <li>RedCell VII*</li> <li>RedCell VIII*</li> <li>RedCell IX*</li> <li>RedCell X*</li> </ol>	Capture-R <sup>®</sup> Select plates

Assay Description	Assay short name	Used reagents	Microplates Used
Antibody identification (0.8 %)	SC_8C_I	1. RedCell 0.8 III* 2. RedCell 0.8 IV* 3. RedCell 0.8 V* 4. RedCell 0.8 VI* 5. RedCell 0.8 VII* 6. RedCell 0.8 VIII* 7. RedCell 0.8 IX* 8. RedCell 0.8 X*	Capture-R <sup>®</sup> Select plates

**\*NOTE:** The reagents described as "RedCell (0.8) III-X" refer to the red cell reagents selected and validated by the customer for use in their laboratory. Immucor will provide 7-digit barcodes for each of the listed RedCell reagents. Those barcodes shall be attached to the selected red cells to enable their use on the NEO Iris<sup>®</sup> and NEO<sup>®</sup> v2.0 instruments. They will not contain any further information on lot numbers and expiry dates. Recording the relevant reagent information will be within the responsibility of the customer. A reagent tracking sheet provided with this Operator Manual Attachment shall facilitate simple tracking of used reagent lots and corresponding expiry dates.

### Limitations, Warnings and Notes

lcon	Description
	<u>Warning:</u> Red cell reagents suitable for the proposed test are required to be selected by the user. Select reagents based on the corresponding red cell Masterlist and further reagent information provided by the reagent's supplier and based on any applicable regulatory requirements. It is the users responsibility to select red cells with the phenotypes matching those specified in local regulations and guidelines applicable to red cell antibody identification techniques.
	<u>Warning:</u> Proper suspension of cellular reagents is essential for unimpeded performance of the assay. For this reason, all cellular reagents, including the customer-chosen red cell reagents (0.8) III – X need to be resuspended by gentle agitation and equipped with a stir ball prior to use. The user needs to make sure that all cellular reagents are placed in the magnetic stirrer positions on the loading bay.
	<u>Warning:</u> The assay SC_8C is designed for the use of <i>RedCell</i> suspensions with a concentration of 2-5% while 0.8% (+/-0.2%) <i>RedCell</i> suspensions can be tested using the SC_8C_I assay. Deviation from the recommended concentrations may lead to false positive, false negative or invalid results. If red blood cells need to be diluted refer to the Instructions for Use of the selected red cells.
	<u>Warning:</u> On board time and reagent expiry of the <i>RedCell</i> reagents cannot be automatically tracked by the system as barcode labels do not contain the relevant information. The operator is required to track lot number, expiry date and first on board time/date for each <i>RedCell</i> vial used on the instrument. A <i>RedCell</i> labeling sheet is provided with this Operator Manual Attachment.
	<u>Warning:</u> Immucor does not provide an on board stability claim for any red cell reagents other than Immucor red cell products released for automated use on the NEO Iris <sup>®</sup> /NEO <sup>®</sup> v2.0. Immucor reagents for automated detection of IgG antibodies on NEO Iris <sup>®</sup> have a 72 hours on board stability claim. Use of Immucor reagent red cells is therefore recommended.

lcon	Description			
	<u>Warning:</u> Reagent fill volume for <i>RedCell</i> reagents is not automatically recognized or tracked in the NEO software. When the <i>RedCell</i> reagent is loaded onto the instrument, the operator is required to manually enter the correct fill volume in the field next to the reagent name. Note that a dead volume of 1 mL needs to be subtracted from the actual fill volume when entering the reagent information, i.e. when the actual vial fill volume is 10 mL, manually type 9 mL in the corresponding data field. To edit the fill volume, click on the relevant loading bay and reagent rack and enter the correct fill volume (see image below):			
	Identifiers Liquid Properties			
	Hack: Prefix + Barcode Settings Volume [m]: Information:			
	1 Reag1 0000421 RedUcel 0.81			
	2 Reag1 0000422 0.000 RedCell 0.8 II			
	3 Reag1 0000423 0.000 RedCell 0.8 III			
	4 Reag1 0000424 0.000 RedCell 0.8 IV			
	<u>Warning:</u> RedCells used on the Identification assays need to be vialed in 10 mL glass vials with the following specifications: 22 x 60 mm, thread 18 mm, flat bottom. Use of vials with differing specifications may lead to abortion of assay runs due to insufficient			
	red cell volume or pipettor crashes.			
	according to the suppliers' specifications to prevent compromising of reagent performance.			
	Warning: Do not pool red cell reagents from different vials or transfer red cell reagents from their original vial to prevent bacterial or chemical contamination.			
	Warning: In a case where both positive control wells obtain reaction scores below the positive cutoff the result for the Capture-R <sup>®</sup> Positive Control is "negative" or in a case where both negative control wells obtain 3+ reactions the result for the Capture-R <sup>®</sup> Negative Control is "positive" due to the definition in the DMS results text. Although the failed control does not display the Ctrl_Fail notice, the run will be invalidated by the system and no results interpretation will be provided.			
	Limitation: It is recommended to use commercially available red cells which have been manufactured specifically for antibody detection and/or identification. Follow the Instructions for Use of the selected red cells.			
	Limitation: It is recommended to place the control strip (first strip) in the first column of the Capture-R <sup>®</sup> Select microplate frame. If the control strip is placed into a different column three negative control results will be displayed on the result reports.			
	<u>Limitation</u> : The 8 cell antibody identification assays are not intended for use with plasma samples originating from Rhlg donors. As samples from this origin are often seen to contain additional antibodies of unknown specifity, false positive results are likely to occur.			
	Limitation: Capture-R Select based Screening Assays do only detect IgG antibodies at a temperature of 37°C.			
	<u>Note:</u> It is recommended to qualify the use of the chosen <i>RedCells</i> by running a sample with a known IgG antibody once per day if the red cell reagents have not been qualified by other QC methods.			



# Description

<u>Note:</u> The RedCell Labeling Sheet that is provided in this Operator Manual Attachment shall help the operator to track relevant *RedCell* reagent data and ensure correct labeling and usage of the red cell products

All further limitations, warnings and notes that relate to the described assays are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings and in Attachment I of the NEO Iris<sup>®</sup> instrument.

# **Cutoffs and Reagent Components**

### List of Cutoffs

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
	Positive/ Negative Control	0	0.0	17.5
		? (not reported)	N/A	N/A
		1+	17.5	50.0
		2+	50.0	72.0
		3+	72.0	90.0
SC_8C		4+	90.0	99.9
SC_8C_I	:_8C_I RedCell (0.8) III/ RedCell III RedCell (0.8) IV/ RedCell IV RedCell (0.8) V/ RedCell V RedCell (0.8) VI/ RedCell VI RedCell VI RedCell VI RedCell VII RedCell VII	0	0.0	20.0
		?	20.0	30.0
		1+	30.0	45.0

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
	RedCell VIII RedCell (0.8) IX/ RedCell IX	2+	45.0	65.0
	RedCell X	3+	65.0	90.0
		4+	90.0	99.9

### Assay Reagent Component Grid

Reagents and Microplates	SC_8C	SC_8C_I
Capture-R <sup>®</sup> Select plate	х	х
corQC <sup>®</sup> EXTEND Standard	х	х
Capture-R <sup>®</sup> Indicator Red Cells	х	х
Capture-R <sup>®</sup> Positive Control Serum (Weak)	х	х
Capture-R <sup>®</sup> Negative Control Serum	х	х
Capture <sup>®</sup> LISS	х	х
RedCell III*	x	
RedCell IV*	x	
RedCell V*	х	
RedCell VI*	х	
RedCell VII*	x	
RedCell VIII*	х	
RedCell IX*	x	

Reagents and Microplates	sc_8c	SC_8C_I
RedCell X*	x	
RedCell 0.8 III*		х
RedCell 0.8 IV*		х
RedCell 0.8 V*		x
RedCell 0.8 VI*		х
RedCell 0.8 VII*		х
RedCell 0.8 VIII*		х
RedCell 0.8 IX*		х
RedCell 0.8 X*		x

**\*NOTE:** The reagents described as "RedCell (0.8) III-X" refer to the red cell reagents selected and validated by the customer for use in their laboratory.

## **Procedural Steps**

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements.

### Information on RedCell reagent labeling


Immucor provides barcode labels to enable the use of customer validated red cell reagents with concentrations of 2-5 % or 0.8 % for antibody identification purposes using the Capture-R Select based 8 Cell Identification assays. The following table lists the assignment of barcodes to the corresponding *RedCell* reagent.

RedCell reagent name	Barcode
RedCell III*	0000413
RedCell IV*	0000414
RedCell V*	0000415
RedCell VI*	0000416
RedCell VII*	0000417
RedCell VIII*	0000418
RedCell IX*	0000419
RedCell X*	0000410
RedCell 0.8 III*	0000423
RedCell 0.8 IV*	0000424
RedCell 0.8 V*	0000425
RedCell 0.8 VI*	0000426
RedCell 0.8 VII*	0000427
RedCell 0.8 VIII*	0000428
RedCell 0.8 IX*	0000429
RedCell 0.8 X*	0000420

To ensure that barcodes are correctly identified by the instrument, they shall only be used on reagent vial that align with the vials specifications described in the "Limitations, Warnings and Notes" section.

The listed barcodes do not automatically provide information on the lot, expiry dates and on board time of the RedCell reagent. For this purpose, a RedCell Labeling Sheet is attached to this Operator Manual Attachment. By filling in this sheet, tracking of lot numbers, expiry dates and further reagent information shall be ensured.

## **Assay Procedural Steps**

This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert.

Assay Button Abbreviation	Brief Synopsis of Assay Procedural Steps
SC_8C SC_8C_J	<ol> <li>Bring all reagents and blood samples to 18–30 °C before testing.</li> <li>Select suitable red cell reagents to be used for antibody identification on the Capture-R Select based Identification assays based on their antigen profile, red cell concentration and suitability to run on the automated NEO Iris<sup>®</sup>/ NEO<sup>®</sup> v2.0 platform.</li> <li>Attach barcode labels provided by Immucor to the chosen red cell reagent vials and track lot number, expiry date and further reagent information using the RedCell Labeling Sheet.</li> <li>Centrifuge the blood samples to separate the plasma from the red blood cells. Remove the caps from the blood sample tubes.</li> <li>Remove the Capture-R<sup>®</sup> Select microplate frame and the desired number of Capture-R<sup>®</sup> Select strips from the pouch.</li> <li>Remove reagent vial caps.</li> <li>Add one striball to each new vial of Capture-R<sup>®</sup> Ready Indicator Red Cells, corQC<sup>®</sup> EXTEND Standard cells and <i>RedCell</i> reagents to be used. Gently agitate each vial to resuspend the red blood cells.</li> <li>Load reagents, microplates, and blood samples onto the NEO Iris following the procedures in Chapter 6 – Instrument Testing Operation.</li> <li>Manually enter the fill volumes of the <i>RedCell</i> reagents by clicking on the lane bay on which the reagents are loaded and manually entering the fill volume minus the vial deadvolume (1 ml) under "Liquid Properties".</li> <li>Assign the SC_8C and/ or SC_8C_l assay to the blood samples, either manually or following the upload worklist procedure.</li> <li>Start the SC_8C and/ or SC_8C_sl assay following the procedures in Chapter 6 – Instrument Testing Operation. The NEO Iris assay, click the Results button on the main menu bar to access the blood sample results.</li> </ol>

## **Results and Interpretations**

#### Introduction

The NEO Iris<sup>®</sup> generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris<sup>®</sup> assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported as Negative, Positive, Equivocal or Invalid. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

#### **Possible Test Well Results**

Possible test well results for assays include:

Indicator	Name	Description
+	Positive	The reaction value is greater than the cutoff value.
-	Negative	The reaction value is less than or equal to the cutoff value.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

#### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations

SC_8C SC_8C_I	+, -, ?, X	<u>Run control interpretations:</u> Positive, Negative, *INV*, Ctrl_Fail <u>Result interpretations:</u> Positive, Negative, No_Int, *INV*
------------------	------------	--

## RedCell Labeling Sheet

Original red cell reagent information		Information for use on NEO Iris				
Red cell reagent product name	Lot number	Cell concentration	Immucor barcode number (enter manually)	Instrument SN	Comments	Operator
	Expiry date	Vial fill volume	Fill volume entered in NEO software? (fill volume - 1ml)	First on board date and time		
			□ (<)			
Red cell reagent product name	Lot number	Cell concentration	Immucor barcode number (enter manually)	Instrument SN	Comments	Operator
	Expiry date	Vial fill volume	Fill volume entered in NEO software? (fill volume - 1ml)	First on board date and time		
			□ (<)			
Red cell reagent product name	Lot number	Cell concentration	Immucor barcode number (enter manually)	Instrument SN	Comments	Operator
	Expiry date	Vial fill volume	Fill volume entered in NEO software? (fill volume - 1ml)	First on board date and time		
			□ (<)			
Red cell reagent product name	Lot number	Cell concentration	Immucor barcode number (enter manually)	Instrument SN	Comments	Operator
	Expiry date	Vial fill volume	Fill volume entered in NEO software? (fill volume - 1ml)	First on board date and time		
			□ (<)			

# **Attachment XX: NEO Iris Operator Manual**

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results, interpretations and interface data used for the AG\_M\_m assay for the detection of M antigen.

Attachment XVI: NEO Iris Operator ManualX	XX-1
Copyrights and DisclaimersX	XX-2
How this Attachment is OrganizedX	XX-4
Description for monoclonal M antigen AssayX	XX-6
Monoclonal M antigen Assay Cutoffs and Reagent Components	X-8
Monoclonal M antigen Assay Procedural StepsXX	X-10
Monoclonal M antigen Assay Results and InterpretationsXX	X-11
Interface Specification InformationXX	X-14
Performance CharacteristicsXX	X-15

# Copyrights and Disclaimers

©2022, Immucor, Inc. NEO Iris<sup>®</sup> is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XVI for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (<u>doccontrol-</u> <u>de@immucor.com</u>).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

# How this Attachment is Organized

#### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

#### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (XVI) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (A-XVI) in parentheses. The first set of characters (NEO Iris\_EU) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

#### Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

#### Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
$\mathbf{\Lambda}$	Warning	Potentially damaging or dangerous
		outcomes if certain critical
		procedural steps are ignored or
		incorrectly executed.
1 îi	Consult instructions for	
	use	

## Description for monoclonal M antigen Assay

#### Introduction

#### Intended Use

Anti-M IgG (Monoclonal) blood group reagent used to detect the M (MNS1) erythorocyte antigen of donors and recipients by direct hemagglutination test on the NEO Iris<sup>®</sup> automated analyzer using the AG\_M\_m assay.

#### List of assays

The following table lists the monoclonal assays for the detection of Anti-M (AG\_M\_m) on the NEO Iris. The AG\_M\_m assay exist in a horizontal and in a vertical configuration.

#### Version 1.00:

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
Monoclonal M	AG_M_m	<ol> <li>Anti-M (Monoclonal)</li> <li>Monoclonal Control</li> </ol>	Untreated microplates
	Blood grouping	Vertical/8 wells per strip	
Monoclonal M	AG_M_m	<ol> <li>Anti-M (Monoclonal)</li> <li>Monoclonal Control</li> </ol>	Untreated microplates

## Limitations, Warnings and Notes

lcon	Description
	<u>Note</u> : The purpose of the Monoclonal Control for AG_M_m assay is to serve as a sample control for Immucor low protein blood grouping reagents. It is expected to indicate those sample related conditions that could lead to spontaneous agglutination with low protein reagents and, therefore, a false positive interpretation for the test. When the Monoclonal Control well yields a positive result, then the instrument would not report the blood type.
	Limitation: Vials of reagents, that have remained continuously on the NEO Iris for 72 hours (3 days) should be removed and replaced with fresh vials. Vials of reagents that are removed from the NEO Iris when not in use and refrigerated can be used up to their expiration dates.
	Limitation: The NEO Iris does not generate an interpretation of mixed-field. Such a mixed-field reaction will be interpreted as positive, negative, or equivocal.
	Limitation for AG_M_m: Samples that exhibit excessive hemolysis or lipemia, or are icteric, should not be tested on the instrument. Hemolytic samples are tested until a concentration of 1065.3 mg/dl without showing erroneous results. Lipemic samples (Tryglyceride) are tested until a concentration of 1674.0 mg/dl without showing erroneous results. Icteric samples (conjugated bilirubin) are tested until a concentration of 47.6 mg/dl without showing erroneous results. Icteric samples (unconjugated bilirubin) are tested until a concentration of 50.0 mg/dl without showing erroneous results

# Monoclonal M Assay Cutoffs and Reagent Components

### List of Cutoffs

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		0	0	25
		? (Query)	25	28
	Anti-M IgG	1+	28	35
	(Monoclonal)	2+	35	50
		3+	50	80
		4+	80	99
AG_M_m	Monoclonal control	0	0	25
		?(Query)	25	25
		1+	25	35
		2+	35	50
		3+	50	80
		4+	80	99

## Assay Reagent Component Grid

Reagents and Microplates	<ul> <li>Untreated Microplates</li> </ul>
For AG_M_m	(barcoded)
	<ul> <li>Specimen Diluent</li> </ul>
	<ul> <li>Monoclonal Control</li> </ul>
	<ul> <li>Anti-M IgG (monoclonal)</li> </ul>

Operator Manual

# Monoclonal M Assay Procedural Steps

## Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button					
Abbreviation	Brief Synopsis of Assay Procedural Steps				
AG M m	1. Bring all reagents and blood samples to 18–30°C before testing.				
	2. Centrifuge the blood samples to separate the plasma from the red blood				
	cells and then remove the caps from those tubes. Process the donor				
	segment blood samples, but do not centrifuge those samples.				
	3. Remove reagent vial caps.				
	4. Load reagents, microplates, and blood samples, onto the NEO Iris following				
	the procedures in Chapter 6 – Instrument Testing Operation.				
	5. Assign the assay (AG_M_m) to the blood samples either manually or				
	following the upload worklist procedure.				
	6. Start the AG_M_m assay following the procedures in Chapter 6 – Instrument				
	Testing Operation. The NEO Iris automatically performs the AG_M_m assay,				
	and records and interprets blood sample results.				
	7. At the completion of the NEO Iris assays, click the Results button on the				
	main menu bar to access the blood samples results.				

# Monoclonal M Assay Results and Interpretations

#### Introduction

The NEO Iris<sup>®</sup> generates a result for each well read by the instrument and an interpretation of the results. The Monoclonal M Assay have predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative, Invalid or Equivocal. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

#### **Possible Test Well Results**

Possible test well results for Monoclonal M antigen Assay include:

Indicator	Name	Description		
+	Positive	The reaction value can be considered as positive.		
-	Negative	The reaction value can be considered as negative.		
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.		
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).		

### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations
AG_M_m	Monoclonal Control: +(1, 2, 3, 4), -, X	Monoclonal Control interpretations:
	M: +(1, 2, 3, 4), -, ? <sup>1</sup> , X	positive, negative, Ctrl failure
		<u>M interpretations:</u> M+, M-, NTD, *INV*

<sup>1</sup> Equivocal results of antibody wells may be edited. There cannot be an equivocal result for Monoclonal Control wells, wherefore negative control results are not editable.

## Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test Interpretation	Description
INV	Invalid
NTD	No Type Determined
Ctrl failure	Control failure
Mono ctrl	Negative control failure
failure	5

# Interface Specification Information

The Monoclonal M assay information for the interface specification, detailed below, is supplemental information to **Attachment II for the NEO Iris Operator Manual (EU).** Refer to **Attachment II for the NEO Iris Operator Manual (EU)** for a full description of the interface specification.

#### **Result Record Field 4 Measurement Values**

SUB-COMPONENT 1				
Assay	Possible Values for Reaction Pattern			
All Assays	-, 1, 2, 3, 4, ?, X			
SUB-COMPONENT 2				
Assay	Interpretation Values			

#### LIS Well Result Identification for Result Field 4 Measurement Values

Assay		Test phase	
	1	Anti-M (Monoclonal)	
AG_M_M	2	Monoclonal Control	

# Performance Characteristics

#### **Specific Performance Characteristics**

#### Sensitivity/Specificity Studies:

Result concordance of Anti-M IgG used on the AG\_M\_m assay on the NEO Iris was compared to results of manual tube test:

The following tables display the result concordance:

#### Anti-M (Monoclonal) Gamma-clone:

Initial Results N=1023		Comparat	tor Reagent		
		Positive	Negative		
	Desitive	805	1*	Positive Percent Agreement	99.9%
Anti-M	Positive			PPA (95% 1-Sided LCI)	99.7%
			217	Negative Percent Agreement	99.5%
	Negative	0		NPA (95% 1-Sided LCI)	97.8%

\* One (1) sample showed discrepancy between M\_PEO and the reference method. The sample remained discrepant after performing genotyping. The root cause of the discrepancy was found to be related to the specimen that was derived from a patient obtaining at least one transfusion before blood drawing which is the propable cause of the inconclusive phenotyping results.

Sensitivity and Specificity were calculated with the initial testing data and displayed in the following table:

<u>Sensitivity</u>	Specificity
<u>100% (805/805)</u>	<u>99.5% (216/217)</u>

#### **Repeatability and Reproducibility Studies:**

Anti-M (Monoclonal) reagent:

Repeatability and Reproducibility were performed with five (5) samples in triplicates, on three (3) intstruments, on five (5) non-consecutive days, two (2) runs per day, over a time period of ten (10) days.. Results demonstrated 100% agreement for all positive and negative test samples

# **Attachment XXI NEO Iris Operator Manual**

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results, interpretations and interface data used for the AG\_N\_m assay for the detection of N antigen.

Attachment XXI NEO Iris Operator Manual	1
Copyrights and Disclaimers	2
How this Attachment is Organized	4
Description for monoclonal N antigen Assay	6
Monoclonal N Assay Cutoffs and Reagent Components	9
Monoclonal N Assay Procedural Steps	11
Monoclonal N Assay Results and Interpretations	12
Interface Specification Information	15
Performance Characteristics	16

# Copyrights and Disclaimers

©2022, Immucor, Inc. NEO Iris<sup>®</sup> is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XVI for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (doccontrol-de@immucor.com).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

# How this Attachment is Organized

#### In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

#### **Notational Conventions**

This attachment uses a page numbering system that includes a prefix of the attachment number (XVI) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (A-XVI) in parentheses. The first set of characters (NEO Iris\_EU) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

#### Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

#### Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to	
$\mathbf{\Lambda}$	Warning	Potentially damaging or dangerous	
		outcomes if certain critical	
		procedural steps are ignored or	
		incorrectly executed.	
1 îi	Consult instructions for		
	use		

## Description for monoclonal N antigen Assay

#### Introduction

#### Intended Use

Anti-N IgM (Monoclonal) blood group reagent used to detect the N (MNS1) erythorocyte antigen of donors and recipients by direct hemagglutination test on the NEO Iris<sup>®</sup> automated analyzer using the AG\_N\_m assay.

#### List of assays

The following table lists the monoclonal assays for the detection of Anti-N (AG\_N\_m) on the NEO Iris. The AG\_N\_m assay exist in a horizontal and in a vertical configuration.

#### Version 1.00:

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
Monoclonal N	AG_N_m	<ol> <li>Anti-N (Monoclonal)</li> <li>Monoclonal Control</li> </ol>	Untreated microplates
	Blood grouping	Vertical/8 wells per strip	
Monoclonal N	AG_N_m	<ol> <li>Anti-N (Monoclonal)</li> <li>Monoclonal Control</li> </ol>	Untreated microplates

## Limitations, Warnings and Notes

lcon	Description
	<u>Note</u> : The purpose of the Monoclonal Control for AG_N_m assay is to serve as a sample control for Immucor low protein blood grouping reagents. It is expected to indicate those sample related conditions that could lead to spontaneous agglutination with low protein reagents and, therefore, a false positive interpretation for the test. When the Monoclonal Control well yields a positive result, then the instrument would not report the blood type.
	Limitation: Vials of reagents, that have remained continuously on the NEO Iris for 72 hours (3 days) should be removed and replaced with fresh vials. Vials of reagents that are removed from the NEO Iris when not in use and refrigerated can be used up to their expiration dates.
	Limitation: The NEO Iris cannot reliably detect hemagglutination reactions that are graded as 1+ or less in test tube methodology. The NEO Iris does not generate an interpretation of mixed-field. Such a mixed-field reaction will be interpreted as positive, negative, or equivocal.
	Limitation: The grading of reactions on the NEO Iris must only be regarded as an approximation when compared to off-line visual grading by laboratory technical staff.

lcon	Description
$\rightarrow$	Limitation for AG_N_m:
<b>_</b> •	Samples that exhibit excessive hemolysis or lipemia, or icteric, should not be tested
	on the instrument. Hemolytic samples are tested until a concentration of 1065.3
	mg/dl without showing erroneous results. Lipemic samples are tested until a
	concentration of 1674.0 mg/dl without showing erroneous results. Icteric samples
	(conjugated bilirubin) are tested until a concentration of 47.6 mg/dl without
	showing erroneous results. Icteric samples (unconjugated bilirubin are tested until a
	concentration of 50.0 mg/dl without showing erroneous results

# Monoclonal N Assay Cutoffs and Reagent Components

## List of Cutoffs

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		0	0	25
	Anti-N IgG (Monoclonal)	? (Query)	25	28
		1+	28	35
		2+	35	50
		3+	50	80
		4+	80	99
AG_N_III		0	0	25
		?(Query)	25	25
	Monoclonal	1+	25	35
	control	2+	35	50
		3+	50	80
		4+	80	99

## Assay Reagent Component Grid

Reagents and Microplates	<ul> <li>Untreated Microplates</li> </ul>
For AG_N_m	(barcoded)
	<ul> <li>Specimen Diluent</li> </ul>
	<ul> <li>Monoclonal Control</li> </ul>
	<ul> <li>Anti-N IgM (monoclonal)</li> </ul>

Operator Manual

# Monoclonal N Assay Procedural Steps

## Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
AG N m	1. Bring all reagents and blood samples to 18–30°C before testing.
	2. Centrifuge the blood samples to separate the plasma from the red blood
	cells and then remove the caps from those tubes. Process the donor
	segment blood samples, but do not centrifuge those samples.
	3. Remove reagent vial caps.
	4. Load reagents, microplates, and blood samples, onto the NEO Iris following
	the procedures in Chapter 6 – Instrument Testing Operation.
	5. Assign the assay (AG_N_m) to the blood samples either manually or
	following the upload worklist procedure.
	6. Start the AG_N_m assay following the procedures in Chapter 6 – Instrument
	Testing Operation. The NEO Iris automatically performs the AG_N_m assay,
	and records and interprets blood sample results.
	7. At the completion of the NEO Iris assays, click the Results button on the
	main menu bar to access the blood samples results.

## Monoclonal N Assay Results and Interpretations

#### Introduction

The NEO Iris<sup>®</sup> generates a result for each well read by the instrument and an interpretation of the results. The Monoclonal N Assay have predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative, Invalid or Equivocal. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

#### **Possible Test Well Results**

Possible test well results for Monoclonal N antigen Assay include:

Indicator	Name	Description
+	Positive	The reaction value can be considered as positive.
-	Negative	The reaction value can be considered as negative.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

#### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations
AG_N_m	Monoclonal Control: +(1, 2, 3, 4), -, X	Monoclonal Control interpretations:
	N: +(1, 2, 3, 4), -, ? <sup>1</sup> , X	positive, negative, Ctrl failure
		<u>N interpretations:</u> N+, N-, NTD, *INV*

<sup>1</sup> Equivocal results of antibody wells may be edited. There cannot be an equivocal result for Monoclonal Control wells, wherefore negative control results are not editable.

### Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test	Description		
Interpretation	Description		
INV	Invalid		
NTD	No Type Determined		
Ctrl failure	Control failure		
Mono ctrl			
failure			
# Interface Specification Information

The Monoclonal N assay information for the interface specification, detailed below, is supplemental information to **Attachment II for the NEO Iris Operator Manual (EU).** Refer to **Attachment II for the NEO Iris Operator Manual (EU)** for a full description of the interface specification.

#### **Result Record Field 4 Measurement Values**

SUB-COMPONENT 1			
Assay	Possible Values for Reaction Pattern		
All Assays	-, 1, 2, 3, 4, ?, X		
SUB-COMPONENT 2			
Assay	Interpretation Values		

#### LIS Well Result Identification for Result Field 4 Measurement Values

Assay	Test phase	
	1	Monoclonal Control
AG_N_m	2	Anti-N (Monoclonal)

# Performance Characteristics

## **Specific Performance Characteristics**

#### Sensitivity/Specificity Studies:

Result concordance of Anti-N IgM used on the AG\_N\_m assay on the NEO Iris was compared to results of manual tube test:

The following tables display the result concordance:

#### Anti-N (Monoclonal) Gamma-clone:

Initial Results		Comparator Reagent				
N=1023	3	Positive Negative				
	Desitive	705	705		Positive Percent Agreement	100%
Anti-N	Positive 705 2*	Ζ*	PPA (95% 1-Sided LCI)	99.6%		
				Negative Percent Agreement	99.4%	
	Negative	0	316	NPA (95% 1-Sided LCI)	98.0%	

\* Two (2) samples were discrepant at initial test. Both were N negative in the tube but positive on the automated method at initial testing. Two (2) re-tests were conducted and have that the initial tube (reference method) test showed a discrepant result. After re-test, the tube test results were N positive in the tube test. Therefore, the discrepancies were resolved.

Sensitivity and Specificity were calculated with the initial testing data and displayed in the following table:

Sensitivity	Specificity
100% (705/705)	99.4% (314/316)

#### **Repeatability and Reproducibility Studies:**

#### Anti-N (Monoclonal) reagent:

Repeatability and Reproducibility were performed with five (5) samples in triplicates, on three (3) intstruments, on five (5) non-consecutive days, two (2) runs per day, over a time period of ten (10) days.. Results demonstrated 100% agreement for all positive and negative test samples

# **Attachment XXII: NEO Iris Operator Manual**

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results, interpretations and interface data used for the AG\_K (H), PHENO12\_2, PHENO16 and PHENO16\_2 assays for the detection of phenotyping status on RhCE and/or Kell.

Attachment XXII: NEO Iris Operator Manual	XXII-1
Copyrights and Disclaimers	XXII-2
How this Attachment is Organized	XXII-4
Description for AG_K, PHENO12_2, PHENO16 and PHENO16_2	XXII-6
Cutoffs and Reagent Components	XXII-8
Procedural Steps	XXII-16
Results and Interpretations	XXII-17
Interface Specification Information	XXII-20
Performance Characteristics	XXII-22

# Copyrights and Disclaimers

©2022, Immucor, Inc. NEO Iris<sup>®</sup> is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XXII for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (<u>doccontrol-de@immucor.com</u>).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

# How this Attachment is Organized

## In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

## Notational Conventions

This attachment uses a page numbering system that includes a prefix of the attachment number (XXII) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (A- XXII) in parentheses. The first set of characters (NEO Iris\_EU) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

## Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

#### Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

lcon	Type of Warning	Related to
$\mathbf{\Lambda}$	Warning	Potentially damaging or dangerous
		outcomes if certain critical
		procedural steps are ignored or
		incorrectly executed.
1 îi	Consult instructions for	
	use	

# Description for AG\_K (H), PHENO12\_2, PHENO16 and PHENO16\_2 assays

## Introduction

#### Intended Use

The described assays are designed to provide phenotyping status on RhCE and/or Kell. Assays are based on direct agglutination technology for the NEO Iris<sup>™</sup> automated analyzer.

#### List of assays

The following table lists the described assays on the NEO Iris. The assays exist in a horizontal or vertical configuration.

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
Rh Blood Group	AG_K	<ol> <li>immuClone Anti-Kell (1)</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates
	Blood grouping	Vertical/8 wells per strip	
Rh Blood Group	PHENO12_2	<ol> <li>immuClone Anti-C(2)</li> <li>immuClone Anti-c(2)</li> <li>immuClone Anti-E(2)</li> <li>immuClone Anti-e(2)</li> <li>immuClone Anti-Kell (1)</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates

#### NEO Iris Operator Manual

	Blood grouping	Horizontal/12 wells per strip	
Rh Blood Group	PHENO16	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Anti-Kell (1)</li> </ol>	Untreated microplates
		6. immuClone Rh-Hr Control	
	Blood grouping	Horizontal/12 wells per strip	
Rh Blood Group	PHENO16_2	<ol> <li>immuClone Anti-C(2)</li> <li>immuClone Anti-c(2)</li> <li>immuClone Anti-E(2)</li> <li>immuClone Anti-e(2)</li> <li>immuClone Anti-Kell (1)</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates

## Limitations, Warnings and Notes

All limitations, warnings and notes that relate to the described assays are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings and in Attachment I of the NEO Iris instrument. No additional limitations that relate to the described assays were identified.

# Cutoffs and Reagent Components

## List of Cutoffs

Assay	Poaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	< =
		0	0	23
		?	23	35
	immuClone	1+(not reported)	N/A	N/A
	Anti-Keii (1)	2+	35	50
		3+	50	80
		4+	80	99
AG_K		0	0	23
	immuClone Rh-Hr Control	? (not reported)	N/A	N/A
		1+	23	35
		2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
		?	23	75
PHENO12_2	ImmuClone Anti-C(2)	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		Grade	>	<=
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-c(2)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone Anti-E(2)	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone Anti-e(2)	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
	immuClone	0	0	23

#### Attachment XXII: NEO Iris Operator Manual

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	< =
	Anti-Kell (1)	?	23	35
		1+(not reported)	N/A	N/A
		2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
		? (not reported)	23	23
	immuClone	1+	23	35
	Kn-Hr Control	2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-C(1)	2+ (not reported)	N/A	N/A
PHENO16		3+	75	80
		4+	80	99
		0	0	23
	ImmuClone	?	23	75
	Anti-c(1)	1+ (not reported)	N/A	N/A

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-E(2)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
	ImmuClone	?	23	75
		1+ (not reported)	N/A	N/A
	Anti-e(2)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	35
	immuClone Anti-Kell (1)	1+(not reported)	N/A	N/A
		2+	35	50
		3+	50	80

#### Attachment XXII: NEO Iris Operator Manual

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		4+	80	99
		0	0	23
		? (not reported)	N/A	N/A
	immuClone	1+	23	35
	KII-HI CONITOI	2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone Anti-C(2)	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
PHENO16_2		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-c(2)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
	ImmuClone	0	0	23
	Anti-E(2)	?	23	75

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		1+ (not	> N/A	<= N/A
		reported) 2+ (not		
		reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-e(2)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	35
	immuClone	1+(not reported)	N/A	N/A
	Anti-Keii (1)	2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
	immuClone	? (not reported)	N/A	N/A
Kn-Hr Control	1+	23	35	
		2+	35	50

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Redection	Grade	>	< =
		3+	50	80
		4+	80	99

# Assay Reagent Component Grid

Reagents and Microplates	AG_K	PHENO12_2	PHENO16	PHENO16_2
Untreated Microplates (barcoded)	х	Х	Х	х
ImmuClone Anti-C(1)			Х	
ImmuClone Anti-C(2)		х		х
ImmuClone Anti-c(1)			х	
ImmuClone Anti-c(2)		х		х
ImmuClone Anti-E(1)			х	
ImmuClone Anti-E(2)		х		х
ImmuClone Anti-e(1)			х	
ImmuClone Anti-e(2)		х		х
ImmuClone Anti-Kell (1)	х	х	х	х
ImmuClone Rh-Hr Control	х	х	х	х

# Procedural Steps

## Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
AG K	1. Bring all reagents and blood samples to 18–30°C before testing.
_	2. Centrifuge the blood samples to separate the plasma from the red blood
PHENO12_2	cells and then remove the caps from those tubes. Process the donor
BUENO44	segment blood samples, but do not centrifuge those samples.
PHENO16	3. Remove reagent vial caps.
PHENO16 2	4. Load reagents, microplates, and blood samples, onto the NEO Iris following
	the procedures in Chapter 6 – Instrument Testing Operation.
	5. Assign the appropriate assay to the blood samples either manually or
	following the upload worklist procedure.
	6. Start the appropriate assay following the procedures in Chapter 6 –
	Instrument Testing Operation. The NEO Iris automatically performs the
	assays, and records and interprets blood sample results.
	7. At the completion of the NEO Iris assay, click the Results button on the main
	menu bar to access the blood samples results.

# Results and Interpretations

#### Introduction

The NEO Iris<sup>®</sup>generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative, Invalid or Equivocal. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

## **Possible Test Well Results**

Possible test well results for assays include:

Indicator	Name	Description
+	Positive	The reaction value can be considered as positive.
-	Negative	The reaction value can be considered as negative.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

## **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations
AG_K	ImmuClone Rh-Hr Control: +(1, 2, 3, 4), -, X	Kell Interpretations: K+, K-,
	ImmuClone Anti-Kell (1): +(1, 2, 3, 4), -, ? <sup>1</sup> , X	NTD / Mixed field ?, NTD, *INV*
PHENO12_2	ImmuClone Rh-Hr Control: +(1, 2, 3, 4), -, X	CE Interpretations: cc ee, CC ee, Cc ee,
PHENO16	Other reagents: +(1, 2, 3, 4), -, ? <sup>1</sup> , X	cc Ee, CC Ee, Cc Ee, cc EE, Cc EE, CC EE,
PHENO16_2		NTD / Mixed field ?, NTD, *INV*
		Kell Interpretations: K+, K-,
		NTD / Mixed field ?, NTD, *INV*

<sup>1</sup> Equivocal results of antibody wells may be edited. There cannot be an equivocal result for ImmuClone Rh-Hr Control wells, wherefore negative control results are not editable.

# Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test	Description
Interpretation	
INV	Invalid
NTD	No Type Determined
NTD / Mixed	No Type Determined ( Miyed field 2
field ?	No Type Determined / Mixed field ?
NTD NTD	No Type Determined No Type Determined

# Interface Specification Information

The AG\_K (H), PHENO12\_2, PHENO16 and PHENO16\_2 assays information for the interface specification, detailed below, is supplemental information to **Attachment II for the NEO Iris Operator Manual (EU).** Refer to **Attachment II for the NEO Iris Operator Manual (EU)** for a full description of the interface specification.

#### **SUB-COMPONENT 1** Assay **Possible Values for Reaction Pattern** -, 1, 2, 3, 4, ?, X AG\_K PHENO12\_2 PHENO16 PHENO16\_2 **SUB-COMPONENT 2 Interpretation Values** Assay **Kell Interpretation** AG\_K K+, K-, NTD, \*INV\*, NTD / Mixed field ? **CE Interpretation Kell Interpretation** PHENO12\_2 CC EE, Cc EE, cc EE, CC Ee, Cc K+, K-, NTD, \*INV\*, NTD / Mixed PHENO16 Ee, cc Ee, CC ee, Cc ee, cc ee field ? PHENO16 2 NTD, \*INV\*, NTD / Mixed field ?

## **Result Record Field 4 Measurement Values**

## LIS Well Result Identification for Result Field 4 Measurement Values

Assay	Test phase		
	1	ImmuClone Anti-Kell (1)	
AG_K	2	ImmuClone Rh-Hr Control	
	1	ImmuClone Anti-C(2)	
	2	ImmuClone Anti-c(2)	
PHENO12_2	3	ImmuClone Anti-E(2)	
	4	ImmuClone Anti-e(2)	
	5	ImmuClone Anti-Kell (1)	

	6	ImmuClone Rh-Hr Control	
	1	ImmuClone Anti-C(1)	
	2	ImmuClone Anti-c(1)	
	3	ImmuClone Anti-E(1)	
PHENOTO	4	ImmuClone Anti-e(1)	
	5	ImmuClone Anti-Kell (1)	
	6	ImmuClone Rh-Hr Control	
	1	ImmuClone Anti-C(2)	
	2	ImmuClone Anti-c(2)	
	3	ImmuClone Anti-E(2)	
	4	ImmuClone Anti-e(2)	
	5	ImmuClone Anti-Kell (1)	
	6	ImmuClone Rh-Hr Control	

# Performance Characteristics

## **Specific Performance Characteristics**

#### Method Comparison Studies:

A method comparison study was performed at one (1) internal site at Immucor Medizinische Diagnostik GmbH. Result concordance of the updated AG\_K (H), PHENO12\_2, PHENO16 and PHENO16\_2 assays (i.e. AGK\_PEO, PH122\_PEO, PH16\_PEO and PH162\_PEO) on the NEO Iris (methods under test) were compared to the validated PHENO12 assay (reference method) on the NEO Iris. To verify the result consistency of the method under test and the reference method the concordance of all test results was calculated by comparing the number of results in agreement with the total number of results. Minitab 20 was used for concordance calculations employing the 1-proportion exact method at the 95% lower bound confidence interval (LCI).

Assay performance represents NEO Iris (software version 3.1 or higher).

Result concordance of AG\_K, PHENO12\_2, PHENO16 and PHENO16\_2 on the NEO Iris/v2.0 were compared to results obtained on PHENO12 on the NEO Iris/v2.0. The following tables display the result concordance.

2x2 Contingency Table		Reference method		Analysis
		Positive	Negative	Number of tests
	Positive	76	0	300 Concordance 95% Lower Confidence interval
Method under test	Negative	0	224	99.2% Positive percent agreement 100% (LCI: 97.0%*) Negative percent agreement 100% (LCI: 99.0%)

#### Data analysis for method under test AGK\_PEO:

\* The PPA (95% 1-sided LCI) was 97.0% due to the low frequency of Kell positive antigens in the Caucausian population (9%). The percent agreement (point estimate) was 100%.

Result Concordance	Result Concordance	
(True positive)	(True negative)	
100% (76/76)	100% (224/224)	

#### Data analysis for method under test PH122\_PEO:

2x2 Contingency Table		Reference method		Analysis
		Positive	Negative	Number of tests
	Positive	76	0	300 Concordance 95% Lower Confidence interval
Method under test	Negative	0	224	99.2% Positive percent agreement 100% (LCI: 97.0%*) Negative percent agreement 100% (LCI: 99.0%)

\* The PPA (95% 1-sided LCI) was 97.0% due to the low frequency of Kell positive antigens in the Caucausian population (9%). The percent agreement (point estimate) was 100%.

Result Concordance	Result Concordance	
(True positive)	(True negative)	
100% (76/76)	100% (224/224)	

#### Data analysis for method under test PH16\_PEO:

2x2 Contingency Table		Reference method		Analysis
		Positive	Negative	Number of tests
	Positive	76	0	300 Concordance 95% Lower Confidence interval
Method under test	Negative	0	224	99.2% Positive percent agreement 100% (LCI: 97.0%*) Negative percent agreement 100% (LCI: 99.0%)

\* The PPA (95% 1-sided LCI) was 97.0% due to the low frequency of Kell positive antigens in the Caucausian population (9%). The percent agreement (point estimate) was 100%.

Result Concordance	Result Concordance	
(True positive)	(True negative)	
100% (76/76)	100% (224/224)	

#### Data analysis for method under test PH162\_PEO:

2x2 Contingency Table		Reference method		Analysis
		Positive	Negative	Number of tests
	Positive	76	0	300 Concordance 95% Lower Confidence interval
Method under test	Negative	0	224	99.2% Positive percent agreement 100% (LCI: 97.0%*) Negative percent agreement 100% (LCI: 99.0%)

\* The PPA (95% 1-sided LCI) was 97.0% due to the low frequency of Kell positive antigens in the Caucausian population (9%). The percent agreement (point estimate) was 100%.

Result Concordance	Result Concordance	
(True positive)	(True negative)	
100% (76/76)	100% (224/224)	

# **Attachment XXIV: NEO Iris Operator Manual**

## In This Attachment

This attachment describes the reagents, cutoff values, basic assay procedural steps, possible well results, interpretations and interface data used for the updated assays using the immuClone (2) Anti-K (Kell) Automated IgM reagent.

Attachment XXIV: NEO Iris Operator Manual	XXIV-1
Copyrights and Disclaimers	XXIV-2
How this Attachment is Organized	XXIV-4
Assay Description	XXIV-6
Cutoffs and Reagent Components	XXIV-11
Procedural Steps	XXIV-44
Results and Interpretations	XXIV-46
Performance Characteristics	XXIV-51

# Copyrights and Disclaimers

©2023, Immucor, Inc. NEO Iris® is a registered trademark of Immucor, Inc. (hereinafter "Immucor").

The contents of this manual are protected by copyright. The Immucor name, logos, related trademarks, and service mark are owned and used in commerce by Immucor and are protected by U.S. and international trademark laws.

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any human or computer language in any form by any means without written permission from the copyright holder.

Unauthorized copying of this publication may not only infringe copyright but also reduce the ability of Immucor to provide accurate and up-to-date information to customers.

No warranties of any nature are extended by this document. You should be very careful to ensure that the use of this information and/or hardware and software material complies with the laws, rules, and regulations of the jurisdictions with respect to which it is used.

All operating instructions must be followed. In no event shall Immucor be held responsible for failures, errors, or other liabilities resulting from a customer's noncompliance with the procedures and precautions outlined in this manual.

The sample screen displays and sample printouts in this Attachment XXIII for NEO Iris Operator Manual are for information and illustration purposes only. Immucor makes no representations or warranties about the accuracy or reliability of the information presented on the screen displays, and this information is not to be used for clinical or maintenance evaluation.

Immucor reserves the right to make changes to the product to improve reliability, function or design, or to discontinue any product at any time without notice or obligation. The material contained in this manual is subject to change without notice. Immucor will not be liable for any consequences resulting from the use of this publication.

Any comments or suggestions regarding this publication should be forwarded to Immucor Medizinische Diagnostik GmbH, Robert-Bosch-Strasse 32, D-63303 Dreieich, Germany (<u>doccontrol-de@immucor.com</u>).

Microsoft, Windows, and the Windows logo are registered trademarks or trademarks of Microsoft Corporation in the United States and other countries. All third-party trademarks, service marks, and trade names are the property of their respective owners and are hereby acknowledged. No responsibility is assumed by Immucor for the use or reliability of software or equipment that is not supplied by Immucor or its affiliated dealers. All warnings and cautions must be reviewed by the Operator prior to using the NEO Iris for the first time.

# How this Attachment is Organized

## In This Section

This section describes the organization of this attachment, including:

- Notational Conventions
- Limitations of Use and Warnings
- Use of Icons

## Notational Conventions

This attachment uses a page numbering system that includes a prefix of the attachment number (XXIV) hyphenated with the page number. Sequential lists that describe step-by-step procedures are included as numbered lists.

The footer of each page contains the current attachment version identified using a nine character hyphenated format combined with a five character attachment alpha-numeric identifier (A-XXIV) in parentheses. The first set of characters (NEO Iris\_EU) identify the instrument. The second set of three characters identifies this document as an attachment to the operator manual (001). The final set of three characters identifies the version of the attachment associated with the operator manual. 100 designates version 1; 200 designates a full version 2 update; and so on.

## Limitations of Use and Warnings

Limitations of use and warnings are located in this attachment, specifically in the body of the text where they are most relevant to the information. An icon draws your attention to limitations of use and warnings.

#### Use of Icons

The following icons appear in this attachment to alert you of warnings or limitations of use.

Icon Type of Warning Related to
---------------------------------

	Warning	Potentially damaging or dangerous outcomes if certain critical procedural steps are ignored or incorrectly executed.
	Limitation	
[]i]	Consult instructions for use	

# Assay Description

## Introduction

#### Intended Use

The described assays are designed to provide basic blood type, Rh status and/or Kell which can include both forward testing (of red blood cell) and reverse testing (of plasma/serum). Assays are based on direct agglutination technology for the NEO Iris<sup>™</sup> and NEO<sup>®</sup> v2.0 automated analyzer.

#### List of assays

The following table lists the described assays on the NEO Iris.

Assay Description	Assay short name	Used reagents	Microplates Used
	Blood grouping	Horizontal/12 wells per strip	
Rh Blood Group	RHFORMEL	<ol> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-C(2)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-c(2)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-E(2)</li> <li>immuClone Anti-e(1)</li> <li>immuClone Anti-e(2)</li> <li>immuClone Anti-Kell (1)</li> <li>ImmuClone (2) Anti-K Automated</li> <li>IgM</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates

Assay Description	Assay short name	Used reagents	Microplates Used	
		1. immuClone Anti-C(1)		
		2. immuClone Anti-C(2)		
		3. immuClone Anti-c(1)		
		4. immuClone Anti-c(2)		
		5. immuClone Anti-E(1)		
Ph Plood		6. immuClone Anti-E(2)	Untroated	
Group	RHFORM_D	7. immuClone Anti-e(1)	microplates	
		8. immuClone Anti-e(2)		
		9. immuClone Anti-Kell (1)		
		10. immuClone (2) Anti-K Automated IgM		
		11. immuClone Anti-D rapid		
		12. immuClone Rh-Hr Control		
	RHFORM_Cw	1. immuClone Anti-C(1)		
		2. immuClone Anti-C(2)		
		3. immuClone Anti-c(1)		
		4. immuClone Anti-c(2)		
		5. immuClone Anti-E(1)		
Ph Blood		6. immuClone Anti-E(2)	Untreated	
Group		7. immuClone Anti-e(1)	microplates	
		8. immuClone Anti-e(2)		
		9. immuClone Anti-Kell (1)		
		10. ImmuClone (2) Anti-K Automated IgM		
		11. immuClone Anti-Cw		
		12. immuClone Rh-Hr Control		

#### Attachment XXIV: NEO Iris Operator Manual

Assay Description	Assay short name	Used reagents	Microplates Used	
Rh Blood Group	PHENO16_3	1. immuClone Anti-C(2)	Untreated microplates	
		2. immuClone Anti-c(2)		
		3. immuClone Anti-E(2)		
		4. immuClone Anti-e(2)		
		5. ImmuClone (2) Anti-K Automated IgM		
		6. immuClone Rh-Hr Control		
Rh Blood Group	PHENO16_4	1. immuClone Anti-C(1)	Untreated microplates	
		2. immuClone Anti-c(1)		
		3. immuClone Anti-E(1)		
		4. immuClone Anti-e(1)		
		5. ImmuClone (2) Anti-K Automated IgM		
		6. immuClone Rh-Hr Control		
Rh Blood Group	KELL	1. immuClone Anti-Kell (1)	Untreated microplates	
		2. ImmuClone (2) Anti-K Automated IgM		
		3. immuClone Rh-Hr Control		
Forward and Reverse ABO Blood Grouping incl. Kell	ABDLONGK	1. NOVACLONE Anti-A	Untreated microplates	
		2. NOVACLONE Anti-B		
		3. NOVACLONE Anti-AB		
		4. A1-Cell		
		5. A2-Cell		
		6. B-Cell		
		7. O-Cell		
		8. Autocontrol		
		9. NOVACLONE Anti-D		
		10. immuClone Anti-D rapid		
		11. ImmuClone (2) Anti-K Automated IgM		
		12. immuClone Rh-Hr Control		
Assay Description	Assay short name	Used reagents	Microplates Used	
-------------------------	------------------	---------------------------------------	------------------	--
		1. immuClone Anti-A		
		2. immuClone Anti-B		
		3. immuClone Anti-AB		
		4. A1-Cell		
Forward and		5. A2-Cell		
Reverse ABO		6. B-Cell		
Blood	ABDLONG2K	7. O-Cell	Untreated	
Grouping incl.		8. Autocontrol	meroplates	
Kell		9. NOVACLONE Anti-D		
		10. immuClone Anti-D rapid		
		11. ImmuClone (2) Anti-K Automated		
		IgM		
		12. immuClone Rh-Hr Control		
		1. ImmuClone Rh-Hr Control		
Forward ABO	ABD6K_I	2. immuClone Anti-A		
Blood		3. immuClone Anti-B	Untreated	
Grouping incl.		4. immuClone Anti-AB	microplates	
Kell		5. immuClone Anti-D rapid		
		6. ImmuClone (2) Anti-K Automated IgM		
		1. NOVACLONE Diluent Control		
Forward ABO		2. NOVACLONE Anti-A		
Blood		3. NOVACLONE Anti-B	Untreated	
Grouping incl.	ABD6K_N	4. NOVACLONE Anti-AB	microplates	
Kell		5. immuClone Anti-D rapid		
		6. ImmuClone (2) Anti-K Automated IgM		
	Blood grouping	Vertical/8 wells per strip		
Rh Blood		1. ImmuClone (2) Anti-K Automated IgM	Untreated	
Group AG_K 2. immuClone		2. immuClone Rh-Hr Control	microplates	

Assay Description	Assay short name	Used reagents	Microplates Used
Rh Blood Group	PHENO	<ol> <li>immuClone Rh-Hr Control</li> <li>immuClone Anti-C(1)</li> <li>immuClone Anti-E(1)</li> <li>immuClone Anti-c(1)</li> <li>immuClone Anti-e(1)</li> <li>ImmuClone (2) Anti-K Automated IgM</li> </ol>	Untreated microplates
Rh Blood Group	PHENO12_3	<ol> <li>immuClone Anti-C(2)</li> <li>immuClone Anti-c(2)</li> <li>immuClone Anti-E(2)</li> <li>immuClone Anti-e(2)</li> <li>ImmuClone (2) Anti-K Automated IgM</li> <li>immuClone Rh-Hr Control</li> </ol>	Untreated microplates

Assay Description	Assay short name	Used reagents	Microplates Used
	QC		
QC	QC	According to customer specific requirements	Untreated microplates

### Limitations, Warnings and Notes

All limitations, warnings and notes that relate to the described assays are defined in the Operator Manual Chapter 12: Limitations of Use and Warnings and in Attachment I of the NEO Iris instrument. No additional limitations that relate to the described assays were identified.

# Cutoffs and Reagent Components

### List of Cutoffs

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-C(1)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
	ImmuClone Anti-C(2)	?	23	75
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
RHFORMEL		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-c(1)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
	ImmuClone Anti-c(2)	?	23	75
	Anu-C(2)	1+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-E(1)	2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
	ImmuClone Anti-E(2)	0	0	23
		?	23	75
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-e(1)	2+ (not reported)	N/A	N/A
		3+	75	80
	4+	80	99	
		0	0	23
	ImmuClone	?	23	75
	/	1+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		2+ (not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	35
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-Kell (1)	2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
	ImmuClone (2) Anti-K Automated	?	23	35
		1+ (not reported)	N/A	N/A
		2+	35	50
	lgM	3+	50	80
		4+	80	99
		0	0	23
		? (not reported)	N/A	N/A
	ImmuClone Rh-Hr Control	1+	23	35
		2+	35	50
		3+	50	76
		4+	76	99

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=

Assay	Poaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-C(1)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone Anti-C(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
RHFORM_D	ImmuClone Anti-c(1)	0	0	23
		?	23	75
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-c(2)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99

Assay	Peaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		0	0.	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(1)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone Anti-E(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
	ImmuClone Anti-e(1)	0	0	23
		?	23	75
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		<u>4</u> +	80	99
		0	0	23
		?	23	35
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-Kell (1)	2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
	ImmuClone (2)	?	23	35
	Anti-K	1+(not reported)	N/A	N/A
	Automated	2+	35	50
	IgM	3+	50	80
		4+	80	99
		0	0	30
		?	30	76
	ImmuClono	1+(not reported)	N/A	N/A
	ImmuClone Anti-D rapid	2+(not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
		0	0	23
	ImmuClone Rh-Hr Control	? (not reported)	N/A	N/A
		1+	23	35
		2+	35	50
		3+	50	76

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		4+	76	99

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-C(1)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
	ImmuClone Anti-C(2)	0	0	23
RHFORM_Cw		?	23	75
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone Anti-c(1)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A

Assay	Peaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-c(2)	2+(not reported)	23     75       N/A     N/A       N/A     N/A       75     80       80     99       0     23       23     75       N/A     N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(1)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	$ \begin{array}{c} \begin{array}{c} \begin{array}{c} \begin{array}{c} \begin{array}{c} \begin{array}{c} \begin{array}{c} \begin{array}{c}$	N/A	
	Anti-E(2)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		0	0	23
	ImmuClone Anti-e(1)	?	23	75
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-e(2)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	35
	ImmuClone Anti-Kell (1)	1+(not reported)	N/A	N/A
		2+	35	50
	3+	50	80	
		4+	80	99
	ImmuClone (2)	0	0	23
	Anti-K	?	23	35

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
	Automated IgM	1+(not reported)	N/A	N/A
		2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone Anti-Cw	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
	ImmuClone	? (not reported)	N/A	N/A
		1+	23	35
		2+	35	50
		3+	50	76
		4+	76	99

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
	ImmuClone	0	0	23
		?	23	75
PHENOTO_5	Anti-C(2)	1+(not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone Anti-c(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		reported)       3+       75       4+       80       0       ?       23       1+(not	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(2)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
	ImmuClone	?	23	75
	Anti-e(2)	1+(not reported)	N/A	N/A

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	35
	ImmuClone (2) Anti-K Automated IaM	1+(not reported)	N/A	N/A
		2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
	ImmuClone Rh-Hr Control	? (not reported)	N/A	N/A
		1+	23	35
		2+	35	50
		3+	50	76
		4+	76	99

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			F	
		0	0	23
PHENO16_4		?	23	75
	ImmuClone Anti-C(1)	1+(not reported) N/A	N/A	N/A
		2+(not reported)	N/A	N/A

Assay	Peaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-c(1)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	80 99 23 75 N/A N/A 80 99 23 75 N/A N/A 80 99 23 75 N/A 80 99 23 75 N/A
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(1)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	23     75       N/A     N/A       N/A     N/A       75     80       80     99       0     23	99
		0	0	23
		?	23	75
	ImmuClone Anti-e(1)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		4+	80	99
		0	0	23
		?	23	35
	ImmuClone (2) Anti-K	1+(not reported)	N/A	N/A
	IgM	2+	35	50
	5	3+	50	80
		4+	80	99
		0	0	23
		? (not reported)	N/A	N/A
	ImmuClone	1+	23	35
		2+	35	50
		3+	50	76
		4+	76	99
Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		0	0	23
		?	23	35
KELL	ImmuClone	1+(not reported)	N/A	N/A
		2+	35	50
		3+	50	80
		4+	80	99
	ImmuClone (2)	0	0	23
	Anti-K	?	23	35

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
	Automated IgM	1+(not reported)	N/A	N/A
		2+	35	50
		3+	50	80
		4+	80	99
		0	0	23
	ImmuClone Rh-Hr Control	? (not reported)	N/A	N/A
		1+	23	35
		2+	35	50
		3+	50	76
		4+	76	99

Assay	Reaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		0	0	30
		? 30 1+ (not reported) N/A	30	58
			N/A	
ABDLONGK	Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported) N/A	N/A	
		4+	58	99
	NOVACLONE	0	0	30
	Anti-B	?	30	76

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation			>	<=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	NOVACLONE Anti-AB	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
		0	0	23
		?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	99
		0	0	23
	A2-Cell	?	23	28
		1+	28	35
		2+	35	50

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Gidde	>	<=
		3+	50	76
		4+	76	99
		0	0	23
		?	23	28
		1+	28	35
	R-Cell	2+	35	50
		3+	50	76
		4+	76	99
		0	0	23
		?	23	28
	O-Cell	1+	28	35
		2+	35	50
		3+	50	76
		4+	76	99
		0	0	23
		?	76 0 23 28 35 50 76 0 23 28 35 50 76 0 23 28 35 50 76 0 23 28 35 50 76 0 23 28 35 50 76 0 23 28 35 50 76 0 23 28 35	28
	Autocontrol	1+	28	35
	Autocontrol	2+	35	50
		3+	50	76
		4+	76	99
		0	0	30
	NOVACLONE	?	30	76
	Anti-D	1+ (not reported)	N/A	N/A

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
		0	0	30
		?	30	76
	ImmuClass	1+ (not reported)	N/A	N/A
	Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
		0	0	23
		?	23	35
	ImmuClone (2) Anti-K	1+ (not reported)	N/A	N/A
	IgM	2+	35	50
		3+	50	80
		4+	80	99
ImmuClone		0	0	30
	ImmuClone	? (not reported)	N/A	N/A
		1+	30	35
		2+	35	50

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		3+	50	76
		4+	76	99

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation		0	>	30
		0	0	50
		?	30	58
	ImmuClone	1+ (not reported)	N/A	N/A
	Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	99
	ImmuClone Anti-B	0	0	30
ABDLONG2K		?	30	76
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
	ImmuClone	0	0	30
	Anti-AB	?	30	76

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
		0	0	23
		?	23	28
		1+	28	35
	AT-Cell	2+	35	50
		3+	50	76
		4+	76	99
		0	0	23
		?	23	28
		1+	28	35
	Az-Cell	2+	35	50
		3+	50	76
		4+	76	99
		0	0	23
	B-Cell	?	23	28
		1+	28	35
		2+	35	50
		3+	50	76
		4+	76	99

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		0	0	23
		?	23	28
		1+	28	35
	O-Cell	2+	35	50
		3+	50	76
		4+	76	99
		0	0	23
		?	23	28
	Autosoutus	1+	28	35
	Autocontrol	2+	35	50
		3+	50	76
	3+	76	99	
		$ \begin{array}{c ccccccccccccccccccccccccccccccccccc$	30	
			76	
		1+ (not reported)	1+ (not reported) N/A	N/A
	NOVACLONE Anti-D	2+ (not reported)	N/A	N/A
	3+ (not reported)	N/A	N/A	
		4+	76	99
	ImmuClone	0	0	30
	Anti-D rapid	?	30	76

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
		0	0	23
		?	23	35
	ImmuClone (2) Anti-K Automated IgM	1+ (not reported)	N/A	N/A
		2+	35	50
		3+	50	80
		4+	80	99
		0	0	30
	ImmuClone Rh-Hr Control	? (not reported)	N/A	N/A
		1+	30	35
		2+	35	50
		3+	50	76
		4+	76	99

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		0	0	30
ABD6K_I	Rh-Hr Control	? (not reported)	N/A	N/A

Assay	Reaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		1+	30	35
		2+	35	50
		3+	50	76
		4+	76	99
		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	ImmuClone Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	99
		0	0	30
		?	30	76
	ImmuClone Anti-B	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
	ImmuClone	0	0	30
	Anti-AB	?	30	76

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
	Anti-D rapid	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
		0	0	23
		?	23	35
ImmuClone (2) Anti-K	1+ (not reported)	N/A	N/A	
	Automated IgM	2+	35	50
		3+	50	80
		4+	80	99

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
ABD6K_N	NOVACLONE	0	0	30

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Diluent Control	? (not reported)	N/A	N/A
		1+	30	35
		2+	35	50
		3+	50	76
		4+	76	99
		0	0	30
		?	30	58
		1+ (not reported)	N/A	N/A
	Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	99
		0	0	30
		?	30	76
		1+ (not reported)	N/A	N/A
NOVACLON Anti-B	Anti-B	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	76	99
	NOVACLONE	0	0	30

Assay	Reaction	Grade	Lower Limit	Upper Limit	
Abbreviation	Redetion	Gidde	>	<=	
	Anti-AB	?	30	76	
		1+ (not reported)	N/A	N/A	
		2+ (not reported)	N/A	N/A	
		3+ (not reported)	N/A	N/A	
		4+	76	99	
		0	0	30	
		?	30	76	
		1+ (not reported)	N/A	N/A	
Anti-D rapid	Anti-D rapid	2+ (not reported)	N/A	N/A	
		3+ (not reported)	N/A	N/A	
		4+	76	99	
		0	0	23	
		?	23	35	
	ImmuClone (2) Anti-K	1+ (not reported)	N/A	N/A	
	IgM	2+	35	50	
	_	3+	50	80	
		4+	80	99	

Assay	Peaction Grade	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=

Assay Abbreviation	Reaction	Grade	Lower Limit	Upper Limit <=
		0	0	23
		?	23	35
ImmuClone (2) Anti-K	ImmuClone (2) Anti-K	1+(not reported)	N/A	N/A
	IgM	2+	35	50
AG_K ImmuClone Rh-Hr Control	3+	50	80	
		4+	80	99
	ImmuClone	0	0	23
		? (not reported)	N/A	N/A
		1+	23	35
	KII-HI CONTO	2+	35	50
		3+	50	80
		4+	80	99

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	< =
	ImmuClone Rh-Hr Control	0	0	23
PHENO		? (not reported)	N/A	N/A
		1+	23	35
		2+	35	50
		3+	50	76
		4+	76	99
	ImmuClone	0	0	23

Assay	Peaction	Grada	Lower Limit	Upper Limit		
Abbreviation	Reaction	Glade	>	<=		
	Anti-C(1)	?	23	75		
		1+(not reported)	N/A	N/A		
		2+(not reported)	N/A	N/A		
		3+	75	80		
		4+	80	99		
		0	0	23		
		?	23	75		
	ImmuClone Anti-E(1)	1+(not reported)	N/A	N/A		
		2+(not reported)	N/A	N/A		
		3+	75	80		
		4+	80	99		
		0	0	23		
		?	23	75		
	ImmuClone	1+(not reported)	N/A	N/A		
	Anti-c(1)	2+(not reported)	N/A	N/A		
		3+	75	80		
		4+	80	99		
	ImmuClone	0	0	23		
	Anti-e(1)	?	23	75		

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
	ImmuClone (2) Anti-K	?	23	35
		1+(not reported)	N/A	N/A
	lgM	2+	35	50
		3+	50	80
		4+	80	99

Assay	Postion	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	< =
		0	0	23
		?	23	75
ImmuClor Anti-C(2 PHENO12_3	ImmuClone Anti-C(2)	1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
	ImmuClone	0	0	23
	Anti-c(2)	?	23	75

Assay	Reaction	Grade	Lower Limit	Upper Limit
Abbreviation	Reaction	Glade	>	<=
		1+(not reported)	N/A	N/A
		2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-E(2)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
		0	0	23
		?	23	75
	ImmuClone	1+(not reported)	N/A	N/A
	Anti-e(2)	2+(not reported)	N/A	N/A
		3+	75	80
		4+	80	99
	ImmuClone (2)	0	0	23
	Anti-K	?	23	35
	Automated IgM		N/A	N/A

Assay	Poaction	Grada	Lower Limit	Upper Limit
Abbreviation	Reaction	Grade	>	<=
		2+	35	50
		3+	50	80
		4+	80	99
	immuClone Rh-Hr Control	0	0	23
		? (not reported)	N/A	N/A
		1+	23	35
		2+	35	50
		3+	50	76
		4+	76	99

### Assay Reagent Component Grid

Reagents and Microplates	RHFORMEL	RHFORM_D	RHFORM_Cw	PHENO16_3	PHENO16_4	KELL	AG_K	PHENO	PHENO12_3
Untreated Microplates (barcoded)	х	х	х	х	х	х	х	х	х
ImmuClone Anti-C(1)	х	х	х		х			х	
ImmuClone Anti-C(2)	х	х	х	х					х
ImmuClone Anti-c(1)	х	х	х		х			х	
ImmuClone Anti-c(2)	х	х	х	х					х
ImmuClone Anti-E(1)	х	х	х		х			х	
ImmuClone Anti-E(2)	х	х	х	х					х
ImmuClone Anti-e(1)	х	х	х		х			х	
ImmuClone Anti-e(2)	х	х	х	х					х
ImmuClone Anti-Kell (1)	х	х	х			х			
ImmuClone (2) Anti-K Automated IgM	х	х	х	х	x	х	х	х	x
ImmuClone Anti-Cw			х						
ImmuClone Anti-D rapid		х							
ImmuClone Rh-Hr Control	х	х	х	х	х	х	х	х	х

Reagents and Microplates	ABDLONGK	<b>ABDLONG2K</b>	ABD6K_I	ABD6K_N
Untreated Microplates (barcoded)	х	х	х	х
Reagent Red Blood cells (Referencells – Group A1)	х	х		
Reagent Red Blood cells (Referencells – Group A2)	х	х		
Reagent Red Blood cells (Referencells – Group B)	х	х		
Reagent Red Blood cells (Referencells – Group O)	х	х		
NOVACLONE Anti-A	х			х
NOVACLONE Anti-B	Х			х
NOVACLONE Anti-AB	х			х
NOVACLONE Anti-D	х	х		
ImmuClone Anti-A		х	х	
ImmuClone Anti-B		х	х	
ImmuClone Anti-AB		х	х	
ImmuClone Anti-D rapid	х	х	х	х
ImmuClone Rh-Hr Control	х	х	х	
NOVACLONE Diluent Control				х
ImmuClone (2) Anti-K Automated IgM	х	х	х	х
Autocontrol	Х	Х		

## Procedural Steps

### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use.

Assay Button	
Abbreviation	Brief Synopsis of Assay Procedural Steps
Brief Synopsis of Assay Procedural Steps	
--	
1. Bring all reagents and blood samples to 18–30°C before testing.	
2. Centrifuge the blood samples to separate the plasma from the red blood	
cells and then remove the caps from those tubes. Process the donor	
segment blood samples, but do not centrifuge those samples.	
3. Remove reagent vial caps.	
4. Load reagents, microplates, and blood samples, onto the NEO Iris following	
the procedures in Chapter 6 – Instrument Testing Operation.	
5. Assign the appropriate assay to the blood samples either manually or	
following the upload worklist procedure.	
6. Start the appropriate assay following the procedures in Chapter 6 –	
Instrument Testing Operation. The NEO Iris automatically performs the	
assays, and records and interprets blood sample results.	
7. At the completion of the NEO Iris assay, click the Results button on the main	
menu bar to access the blood samples results.	

# Results and Interpretations

#### Introduction

The NEO Iris<sup>®</sup> generates a result for each well read by the instrument and an interpretation of the results. Each NEO Iris assay has predefined interpretations for test well results.

A test well result is the reaction result for a given test well. Test well results are reported Positive, Negative, Invalid or Equivocal. These results are determined by comparing the well reaction value to assay-specific cutoff values. Assay cutoff values are listed in this attachment.

#### **Possible Test Well Results**

Possible test well results for assays include:

Indicator	Name	Description
+	Positive	The reaction value can be considered as positive.
-	Negative	The reaction value can be considered as negative.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

#### Results and Interpretations

#### Manual

The described assays are designed to be read automatically by the instrument, and even though you cannot edit reactions or results, you can still view the images of the assay reactions by using the instrument software. Reaction features are described in the associated reagent package inserts.

#### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well	Possible Test Interpretations
	Results	
RHFORMEL	ImmuClone Rh-Hr Control: +(1, 2, 3, 4), -, X Other reagents: +(1, 2, 3, 4), -, ? <sup>1</sup> , X	<u>CE Interpretations:</u> cc ee, Cc EE, CC EE, NTD / Mixed field ?, NTD, *INV* <u>Kell Interpretations:</u> K+, K-, NTD / Mixed field ?, NTD, *INV*

RHFORM_D	+, -, ?, X	CDE Interpretations:	ccD.ee, CCD.ee, CcD.ee, ccD.Ee, CCD.Ee, CcD.Ee,
			ccD.EE, CcD.EE, CCD.EE, ccddee, CCddee,
			Ccddee, ccddEe, CCddEe, CcddEe, ccddEE,
			CcddEE, CCddEE, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*
RHFORM_Cw	+, -, ?, X	CE and Cw Interpretati	ions: cc ee Cw+, CC ee Cw+, Cc ee Cw+, cc
			Ee Cw+, CC Ee Cw+, Cc Ee Cw+, cc EE
			Cw+, Cc EE Cw+, CC EE Cw+, cc ee Cw-
			CC ee Cw-, Cc ee Cw-, ccEe Cw-, CC Ee
			Cw-, Cc Ee Cw-, cc EE Cw-, CcEE Cw-,
			CC EE Cw-, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*
RHFORM_Cw	+, -, ?, X	CE Interpretations:	cc ee, CC ee, Cc ee, cc Ee Cw+, CC Ee,
*Alternative			Cc Ee, cc EE, Cc EE, CC EE,
Aurora File			
			NTD / Mixed field ?, *INV*
		Cw Interpretations:	Cw+, Cw-, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*
PHENO16_3	+, -, ?, X	CE Interpretations:	CC EE,
PHENO16_4			Cc ee, cc ee, NTD / Mixed field ?, *INV*
PHENO12_3		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*
PHENO	+, -, ?, X	CE Interpretations:	CC EE,
			Cc ee, cc ee, IND / DP, *INV*
		Kell Interpretations:	Kell +, Kell -, IND / DP, *INV*
KELL	+, -, ?, X	K+, K-, NTD / Mixed fi	eld ?, *INV*
AG_K			

Attachment XXIV: NEO Iris Operator

Manual
iviuliuui

ABDLONGK	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD/Mixed field ?, *INV*
ABDLONG2K		Rh Interpretations:	RH+, RH-, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*
ABD6K_I	+, -, ?, X	ABO Interpretations:	A, B, AB, O, NTD/Mixed field ?, *INV*
ABD6K_N		Rh Interpretations:	RH+, RH-, NTD / Mixed field ?, *INV*
		Kell Interpretations:	K+, K-, NTD / Mixed field ?, *INV*

<sup>1</sup> Equivocal results of antibody wells may be edited. There cannot be an equivocal result for ImmuClone Rh-Hr Control wells, wherefore negative control results are not editable.

\*Alternative result interpretation depending on installed aurora file version

#### Key for Table of Well Results and Test Interpretations by Assay

The two (2) tables below describe the possible well test results per the Result Interpretations section of this attachment.

Well Result	Description
+	Positive
-	Negative
?	Equivocal
Х	Invalid

Test	Description
Interpretation	Description
INV	Invalid
NTD	No Type Determined
NTD / Mixed	No Type Determined ( Mixed field 2
field ?	No Type Determined / Mixed field ?
IND / DP	No Type Determined / Mixed field ?

NTD NTD No Type Determined No Type Determined

# Performance Characteristics

#### **Specific Performance Characteristics**

#### Method Comparison Studies:

A method comparison study was conducted internally at Immucor Inc. to evaluate the performance of the immuClone(2) Anti-K (Kell) Automated IgM reagent.

Result concordance of the updated assays was determined by verifying that the results generated on the NEO Iris and Galileo NEO System (software version 3.1 or higher) with assays using the immuClone(2) Anti-K (Kell) Automated IgM are equivalent to the results generated for the same samples tested in reference assays using immuClone(1) Anti-K (Kell). Both vertical and horizontal assay configurations were tested for equivalency to ensure representation of all assays that are configured to use the immuClone(2) Anti-K Automated IgM product.

The following acceptance criterion was met for the immuClone (2) Anti-K (Kell) Automated IgM reagent within all assays under test (AG\_K, PHENO12\_3, ABDLONG2K, Kell, PHENO16\_4, RHFORMEL):

• Result interpretations for phenotyping reactions by the System shall be at least 99% (PE) overall concordant, 99% (PE) PPA concordant, and 99% (PE) NPA concordant to the expected result of the test well.

Reagent Concordance Summary for immuClone (2) Anti-K						
Overall Concordance		РРА		NPA		
Assay	PE	Test Count	PE	Test Count	PE	Test Count
AG_K	100%	301/301	100%	109/109	100%	192/192
PHENO12_3	100%	303/303	100%	109/109	100%	194/194
ABDLONG2K	100%	300/300	100%	107/107	100%	193/193
Kell	100%	309/309	100%	108/108	100%	201/201
PHENO16_4	100%	309/309	100%	108/108	100%	201/201
RHFORMEL	100%	309/309	100%	108/108	100%	201/201

Sensitivity and specificity were calculated from all of the test wells generated with the immuClone(2) Anti-K product compared to the true result.

Sensitivity	Specificity
100% (649/649)	100% (1182/1182)

#### **Precision Studies**

#### **Reproducibility and Repeatability Studies:**

This reproducicility and repeatability study was conducted internally at Immucor Inc. to evaluate performance of the immuClone(2) Anti-K (Kell) Automated IgM reagent using the vertical AG\_K assay on the NEO Iris and Galileo NEO System (software version 3.1 or higher). Three (3) samples were tested in quadruplicate (12 total replicates) on four (4) consecutive days and one (1) nonconsecutive day spanning a total of six (6) days.

Twelve (12) replicates/tubes were run with the AG\_K assay two (2) times per day with at least two (2) hours separating the two (2) runs on three (3) NEO Iris / Galileo NEO Systems.

All expected-positive samples tested in replicates generated positive results, and all expectednegative samples tested in replicates generated negative results with the AG\_K assay using the immuClone(2) Anti-K (Kell) Automated IgM reagent. Results generated met all acceptance criteria and demonstrated 100% agreement:

Agreement Analysis	Concordance
Day	100% (PE)
System	100% (PE)
Run (by system)	100% (PE)
Sample	100% (PE)
РРА	100% (PE)

Agreement Analysis	Concordance
NPA	100% (PE)
Overall concordance	100% (PE)

# **NEO Iris<sup>®</sup> ASSAY GUIDE:**

# **Open Channel Anti-D Group Assays**

# In This Chapter:

NEO Iris® ASSAY GUIDE:	1
Open Channel Anti-D Group Assays	1
About this Assay Guide	2
Intended Use	2
List of Assays and Reagents for the NEO ${\sf Iris}^{\circledast}$	2
Assay Procedural Steps	5
Before you begin	5
Test Results and Interpretations	7
Possible Test Well Results	7
Result Interpretations	7
Assay Cutoffs for NEO $Iris^{\circledast}$	g
Reagent Barcodes	10
Limitations and Warnings	12
Example Labelling Sheet for Open Channel Assays	13

# About this Assay Guide

This assay guide is supplementary instructions for use to the NEO Iris<sup>®</sup> Operation Manual. The guide provides specific information about the intended use, operation instructions, assay parameters, limitations and warnings that are unique to the Open Channel Anti-D Assays. Users must still comply with the general instructions provided by the NEO Iris<sup>®</sup> Operator manuals and corresponding reagent package inserts when using this assay, unless others indicated within this guide.

# Intended Use

The Open Channel Anti-D Assays on the NEO Iris<sup>®</sup> fully automates standard direct microplate-based methods for ABO Grouping and Rh(D) typings. The assays are intended to provide the customer the ability to choose alternative Anti-D reagents, potentially beyond those provided by Immucor. The assays employ other Immucor reagents. By leaving a single open channel explicitly for customer selection of an Anti-D reagent, customers may choose a reagent that allows detection of specific Rh(D) categories of special interest (for example, D Category VI). Therefore, depending on the specifics of which Anti-D reagent is selected by the customer, the assay(s) may only be used for donor and/or patient testing as indicated and validated.

The NEO Iris<sup>®</sup> OCABDCHK2 and OCABD6\_I assays are intended to be used for confirmation of ABO and RhD antigens on human erythrocytes from previously characterized samples (patient and/or donor in origin), as long as the end user selects an Anti-D reagent manufactured for and verified by the customer upon the intended population from which the samples arise.

The NEO Iris<sup>®</sup> OCABDLNG2 assay is intended to be used for characterization of ABO and RhD antigens on human erythrocytes of samples (patient and/or donor in origin), as long as the end user selects an Anti-D reagent manufactured for and verified by the customer upon the intended population from which the samples arise.

The ABO portion of the Open Channel assays are intended for use only with the reagents described in **List of Assays section.** The reagent described as OC Anti-D refers to the Anti-D reagent that the customer selects and fully validates its use for each and every customer selected Anti-D reagent prior to any use in their laboratory.

# List of Assays and Reagents for the NEO Iris®

Assay Description	Assay Short Name	Used Reagents	Microplates used
	Open channel	Horizontal/4 wells per strip	
	(Agglutination)		
Forward ABO	OCABDCHK2	1. immuClone Rh-Hr	Untreated microplates
Blood Grouping		Control	
		2. immuClone Anti-A	
		3. immuClone Anti-B	
		4. OC Anti-D	
	Open channel	Horizontal/6 wells per strip	
	(Agglutination)		
Forward ABO	OCABD6_I	1. immuClone Rh-Hr	Untreated microplates
Blood Grouping		Control	
		2. immuClone Anti-A	
		3. immuClone Anti-B	
		4. immuClone Anti-AB	
		5. OC Anti-D	
		6. Novaclone Anti-D	
	Open channel	Horizontal/12 wells per strip	
	(Agglutination)		
Forward and	OCABDLNG2	1. immuClone Rh-Hr	Untreated microplates
reverse ABO		Control	
Blood Grouping		2. immuClone Anti-A	
		3. immuClone Anti-B	
		4. immuClone Anti-AB	
		5. OC Anti-D	
		6. NOVACLONE Anti-D	
		7. Referencells A1	
		8. Referencells A2	
		9. Referencells B	
		10. Referencells O	



#### Warning:

The Open Channel Anti-D Assay is a direct hemagglutination method. When selecting an Anti-D for use in the open channel, a low protein anti-D reagent capable of direct agglutination should be selected. Reagents that are not capable of direct agglutination may give unexpected negative results. High protein reagents may lead to false positive reactions that are not detected by the Rh-Hr control.

The degree to which reagent Anti-D clones can detect D antigen variants and quantitative weak D antigens in a direct test varies greatly. The open channel Anti-D assays are direct test methods. Results of an indirect test should be evaluated before assigning an RhD negative status.

# Assay Procedural Steps

#### Before you begin



You must prepare all of the necessary reagents and samples for each assay according to the detailed reagent package insert requirements. This brief synopsis of assay procedural steps relating to sample and reagent preparation is intended as summarized steps only and is not intended as a substitute for the detailed package insert. Package inserts are also the source of information for limitations of the reagents in use. To confirm the correct reactivity of the reagents, it is recommended that the reagents are to be tested each day of use by running the Quality Control (QC). For QC frequency, minimum requirements refer to national guidelines.

Assay Button				
Abbreviation	Brief Synopsis of Assay Procedural Steps			
OCABDCHK2	Customers must validate these assays and document their validation prior to			
OCABD6_I	employing them for clinical use!			
OCABDLNG2				
	1. Bring all reagents, and blood samples to 18–30°C before testing			
	2. Centrifuge the blood sample collection tubes to separate the plasma from the			
	red blood cells and then remove the caps from those tubes.			
	3. Select the necessary number of untreated microplates for use.			
	4. For OC Anti-D reagent only: Before first use of a vial, stick a new barcode			
	provided by Immucor to the vial (NOTE: Only vials that fit onto the NEO Iris			
	instrument are allowed to be used. Contact your local Application Specialist for			
	further information.) Document the labelling on the open channel labelling			
	sheet.			
	5. Remove reagent vial caps.			
	6. If needed, add one stirball to each new vial of Referencells®. Gently agitate			
	each vial to resuspend the red blood cells.			
	7. Load reagents, microplates, and blood samples onto the NEO Iris following the			
	procedures in Chapter 6 – Instrument Testing Operation.			
	8. Assign the corresponding assay to the blood samples, either manually or			
	following the upload worklist procedure.			
	9. Start the corresponding assay following the procedures in Chapter 6 –			
	Instrument Testing Operation. The NEO Iris automatically performs the assay,			
	and records and interprets the blood sample results.			
	10. At the completion of the assay, click the Results button on the main menu bar			
	to access the blood sample results.			

# Test Results and Interpretations

#### Possible Test Well Results

Possible test well results for all assays include:

Indicator	Name	Description
+	Positive	The reaction value is greater than the positive cutoff value.
-	Negative	The reaction value is less than or equal to the negative cutoff value.
?	Equivocal	The Equivocal symbol indicates that the reaction well cannot be definitively considered negative or positive. The reaction value is greater than the negative cutoff value or equal to or less than the positive cutoff value.
X	Invalid	The Invalid symbol indicates an error status for a given well. An Invalid test well result is generated if the instrument detects a processing error or a process control parameter out of range (e.g. liquid level detection error, clot detection, incubator temperature out of range, etc.).

#### **Result Interpretations**

The interpretation of the test well results is based on the reaction or reaction pattern of individual test well results, applicable control well results, and NEO Iris<sup>®</sup> process control monitoring. Possible interpretations generated by each assay are listed in the table below.

Assay	Possible Well Results	Possible Test Interpretations			
OCABDCHK2		ABO Interpretations:	A, B, AB, O, NTD / Mixed Field?, *INV*		
OCABD6_I	+, - , ?, X	Rh Interpretations:	Positive, Negative, NTD / Mixed Field?,		
OCABDLNG2		*INV*.			

<u>Note</u>: Assays that include the use of two different Anti-D reagents and yield a positive result with one Anti-D reagent, but a negative result with the other Anti-D reagent for a given sample, will generate an NTD result related to the Rh (D) testing. This can be due to, but is not limited to, reduced D antigen epitope representation on the red blood cell membrane, or the presence of a D variant. Clinical determinations are the purview of the lab director, physicians involved in the care of the patient(s) and should fully understand the work up, including all limitations to any selected tests, that may be required as confirmation, before any result is provided or clinical decisions made regarding the result.

# Assay Cutoffs for NEO Iris®

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
OCABDCHK2		0	0	30
		?	30	58
	immu Clana Anti A	1+ (not reported)	N/A	N/A
	Immucione Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
		0	0	30
	immuClone Anti-B	?	30	76
	OC Anti-D	1+(not reported)	N/A	N/A
	immuClone Rh-Hr	2+(not reported)	N/A	N/A
	Control	3+(not reported)	N/A	N/A
		4+	76	100

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
OCABD6_I		0	0	30
		?	30	58
	immuClana Anti A	1+ (not reported)	N/A	N/A
		2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
	immuClone Anti-B	0	0	30
	immuClone Anti-A,B	?	30	76
	OC Anti-D,	1+(not reported)	N/A	N/A
	Novaclone Anti-D	2+(not reported)	N/A	N/A
	immuClone Rh-Hr	3+(not reported)	N/A	N/A
	Control	4+	76	100

Assay Abbreviation	Reaction	Grade	Lower Limit >	Upper Limit <=
OCABDLNG2		0	0	30
		?	30	58
	immulana Anti A	1+ (not reported)	N/A	N/A
	Immucione Anti-A	2+ (not reported)	N/A	N/A
		3+ (not reported)	N/A	N/A
		4+	58	100
	immuClone Anti-B	0	0	30
	immuClone Anti-A,B	?	30	76
	OC Anti-D,	1+ (not reported)	N/A	N/A
	Novaclone Anti-D	2+ (not reported)	N/A	N/A
	immuClone Rh-Hr	3+ (not reported)	N/A	N/A
	Control	4+	76	100
		0	0	23
	Referencells A1, A2, B,	?	23	28
	0	1+	28	35
	Autocontrol	2+	35	50
		3+	50	76
		4+	76	100

# Reagent Barcodes

These NEO Iris<sup>®</sup> assays use barcodes to identify reagents and samples in the loading modules. Using barcode technology increases the number of steps that can be automated, thus decreasing handling errors. Immucor manufactured reagents are provided with a barcode.

#### Immucor Standard Reagent Barcodes

Immucor manufactured reagents are provided with a barcode that is encoded with the following information:

- Reagent ID The type of reagent
- Expiration date The last date that the reagent may be used
- Lot number Batch identifier
- Serial number Unique identifier for each vial

The serial number identifies each reagent vial as unique and allows the system to track consumption of the reagent during use on the system. This provides a process control to monitor fill levels of the vial. The system assumes vials with 10mL reagent are being used.

#### **Customer Selected OC Anti-D Reagent**

Since the OC Anti-D is selected by the customer and may not bear a barcode, or a barcode that is formatted to meet the NEO Iris<sup>®</sup> system requirements, Immucor has available a secondary barcode label that can be applied to the customer validated OC Anti-D Reagent. The barcode provided includes the barcode elements described above with the following key limitations:

- **Reagent ID** an ID character corresponding to the OC Anti-D
- **Expiration date** A series of barcodes will possess an expiration date corresponding to the last day of the indicated month and year.
- Lot number This will be a sequential number, but will not be the exact lot number applied by the reagent manufacturer. Immucor has provided a log sheet as a part of this assay guide so that end users can maintain documented between the NEO recognized lot and the manufacturers lot.
- Serial number Unique identifier for each vial

# Limitations and Warnings

The specific performance characteristics are described in the package inserts of each reagent. The expected results are specific to the reagent or test wells in use.

The limitations and warnings published in the NEO Iris<sup>®</sup> Operation Manual must be observed when using these assays.



#### Warning:

As the Open Channel Anti-D is designed to allow customer selection of an appropriate Anti-D reagent for their use, Immucor has not fully evaluated any specific anti-D performance characteristics for reagents not labeled and sold by Immucor, and as such, cannot provide any assurances as to functionality of the reagent or the combined reagent assay pairing.

The end user must validate open channel assays for their use based upon the customer's intended purpose, Immucor cannot be responsible for the choice of reagent. The open channel assays are to be used in combination with the open channel label provided by Immucor as described. Labelling of Anti-D reagents shall be documented on the open channel labelling sheet (see below) or in laboratory determined methodology. Only vials that fit onto the NEO Iris<sup>®</sup> may be used (contact your local application specialist for hardware support) and Immucor does not provide any guidance other than a 10 ml vial compatible with the instrument rack system.

As with all reagents, it is necessary to avoid contaminating the reagents during use. Signs of potentially unusable reagents include but are not limited to: marked turbidity, precipitate, fibrin gel or particles present. Do not use contaminated reagents. Do not use leaking vials. Vsiually inspect reagent before use. Discard reagent if there is any concern.

# Example Labelling Sheet for Open Channel Assays

The specific performance characteristics are described in the package inserts of each reagent or test well. The expected results are specific to the reagent or test wells in use.

Reagent	Lot no	Evning data	Labelling date	NEO Iric® barcada	Date/Signature	Barcode expiry matches
name	Lot no.	Expiry date	Labelling date	NEO INS <sup>®</sup> barcode		reagent expiry? (Y/N)

# Index

#### In This Index

This index contains an alphabetical listing of subject matter referenced in the NEO Iris Operator Manual with associated page numbers. This enables you to quickly locate specific information about the NEO within this manual.

### 1

14-lane and 5-lane bays barcodes not interpreted, 11-60 color scheme, 1-5 safety features, 2-19 technical data, C-6 troubleshooting, 11-60
14-lane bay dialog, 3-75 Identifiers area, 3-77 parts, 3-76 Rack area, 3-76 Test Selection area, 3-77

#### 5

5-lane bay dialog, 3-79 parts, 3-79 Rack area, 3-80 Reagent IDs area, 3-80 Reagent Properties area, 3-80

### A

ABO Anti-A,B reagent quality control (QC) completing, 10-21 ABO reagent quality control (QC) completing, 10-20 accessing plate based reports, 8-25 sample based reports, 8-27 adding users, 4-3 additional sample results

viewing, 7-15 ambient light excessive, A-5 approving results, 7-17 archive access specific tab, 4-21 configuration, 4-18 creating, 10-32 general options tab, 4-20 schedule next archive tab, 4-20 SQL server specific tab, 4-21 viewing, 7-20 viewing logs, 4-22 assay selections downloading from LIS, 6-14 assays cutoffs, I-1, I-24, II-1 descriptions, I-8 manually selecting, 6-9 procedural steps, I-51 reagent component grid, I-48 reagents, I-1, II-1 requesting additional, 6-36 assigning passwords, 4-2 user access rights, 4-2

### B

barcode settings, 6-5 barcodes, 1-6, 10 microplates, 1-8 micro-well strip, 11 reagent, 10 reagents, 1-8 sample, 1-9 scanning unread, 6-7 symbologies, 1-6 biohazards potential, A-4

## С

cabinet description, 2-5 exterior, 2-8 interior, 2-5 safety features, 2-6 technical data, C-4 camera lamp changing, 10-117 camera reader description, 2-36 door not opening, 11-58 safety features, 2-36 technical data, C-11 troubleshooting, 11-58 centrifugation stalling, 11-39 centrifuge description, 2-31 safety features, 2-34 technical data, C-10 troubleshooting, 11-32 changing camera lamp, 10-117 checking pipettor reference, 10-17

reader performance, 10-29 residual volume, 10-69 cleaning common waste container, 10-43 instrument, 10-15 Liquid Level Sensor and Trough, 10-88 Liquid Overflow Detection Mat, 10-98 LLS, 10-90 pipettor wash towers, 10-40 reader mirror and light diffuser, ... clearing test data, 10-47 clot detection recovery, 11-18 color code LED, 1-3 color scheme 14-lane and 5-lane bays, 1-5 loading tower, 1-3 **Common Buttons** description, 3-5 common waste container cleaning, 10-43 draining, 6-40 emptying, 10-12 completing ABO Anti-A,B reagent quality control (QC), 10-21 ABO reagent quality control (QC), 10-20 kell phenotyping quality control (QC), 10-22 QC3\_Cell quality control (QC), 10-22 sample loading process, 6-16 configuration archive, 4-18 software, A-13

**Confirmation Dialogs** dimensions, C-2 description, 3-5 general, C-2 incubation time, C-9 connections incubator, C-9 instrument, A-9 microplate internal barcode scanner, C-5 descriptions, A-9 personal computer, C-3 making, A-9 pipetting system, C-7 PC, A-10 pipettors, C-7 descriptions, A-10 plate loading tower, C-4 technical data, C-2 power requirements, C-2 consistent color code, 1-3 pump and liquid unit, C-4 continual access, 1-14 pump rack, C-7 continuous loading reagent loading bay, C-6 during operation, 6-36 sample/reagenet barcode scanner, C-6 software, C-3 controls transport system, C-4 loading, 6-21 washer, C-10 replenishing, 6-39 weight, C-2 creating decontaminating archives, 10-32 tubings, 10-49 crossmatch default access rights worklists, II-32 users, 4-14 **Current Report** deleting body, 8-15 historic plate entries, 10-23, 10-48 current result reports users, 4-12 overview, 8-2 description cutoffs 14-lane and 5-lane bays, 2-18 assay, I-24 cabinet, 2-5 camera reader, 2-36 П centrifuge, 2-31 data hood, 2-9 hardware technical, C-1 incubator, 2-26 14-lane and 5-lane bays, C-6 loading tower, 2-11 cabinet, C-4 personal computer, 2-3 camera reader, C-11 pipetting system, 2-20 centrifuge, C-10 plate carrier, 2-17

connections, C-2

transport system, 2-14 washer, 2-28 descriptions assays, 1-8 dimensions technical data, C-2 disposal liquid waste, 2-38 solid waste, 2-38 downloading assay selections from LIS, 6-14 requests from LIS, 6-14 draining common waste container, 6-40 dynamic schedule, 1-13

# E

editing users, 4-10 emptying common waste container, 10-12 system liquid container, 10-63 entries deleting historic plate, 10-23, 10-48 environmental conditions, A-5 ambient light excessive, A-5 humidity excessive, A-5 error codes troubleshooting, 11-4 exiting results, 7-19 test details, 7-16

expected results, 1-14 exporting results, 7-18 extended shutdown NEO Iris, 9-9

# F

failures clot detection, 11-18 pipettor self check, 11-14 software, 11-10 filling system liquid, 10-8 flushing system liquid, 10-62

### G

Galileo Echo intended use, 2

# Η

hardware techincal data general, C-2 technical data, C-1 14-lane and 5-lane bays, C-6 cabinet, C-4 camera reader, C-11 centrifuge, C-10 connections, C-2 dimensions, C-2 incubation time, C-9 incubator, C-9 microplate internal barcode scanner, C-5 personal computer, C-3

pipetting system, C-7 pipettors, C-7 plate loading tower, C-4 power requirements, C-2 pump and liquid unit, C-4 pump rack, C-7 reagent loading bay, C-6 sample/reagent barcode scanner, C-6 software, C-3 transport system, C-4 washer, C-10 weight, C-2 understanding, 2-1 Help accessing, 3-52 blud\_direct, 3-55 options, 3-56 dialog, 3-52 Extended Help, 3-53 Version Information, 3-54 historic plate entries deleting, 10-23, 10-48 hood description, 2-9 safety features, 2-10 humidity excessive, A-5

### /

icons main menu bar, 3-8 test results, 7-7 incubation time technical data, C-9 incubator description, 2-26

dialog, 3-61 door not closing, 11-45 door not opening, 11-46 safety features, 2-27 technical data, C-9 troubleshooting, 11-44 wrong door opening, 11-48 initializing the instrument, 3-11 inspecting syringes, 10-45 installation validating, A-15 installing software, A-12 instrument cleaning, 10-15 connections, A-9 descriptions, A-9 Instrument Settings dialog, 3-46 Barcode Settings Configuration, 3-47 Creating Assay Profiles, 3-50 instrument start-up logging in, 5-3 NEO Iris, 5-1 Starting UP, 5-2 intended use, 1-2 Galileo Echo, 2 internal voltage, A-4 introduction to the NEO Iris, 1-1

#### K

kell phenotyping quality control (QC) completing, 10-22

Keyboard description, 3-6

#### L

layout screen, 3-2 limitations of use, 12-2 Liquid Level Sensor and Trough cleaning, 10-88 Liquid Overflow Detection Mat cleaning, 10-98 liquid waste disposal, 2-38 removing, 6-39 LLS cleaning, 10-90 reinserting, 10-94 removing, 10-90 loading additional samples, 6-36 continuous during operation, 6-36 plates, 6-27 reagents and controls, 6-21 samples, 6-3, I-5 barcode settings, 6-5 procedure, 6-3 loading bays 14-lane and 5-lane, 2-18 loading tower color scheme, 1-3 description, 2-11 safety features, 2-13 logging in, 3-9, 5-3 logging out, 9-2

#### Μ

machine monitor, 3-59 description, 3-3 parts, 3-59 main menu bar, 3-8 description, 3-3 icons, 3-8 main screen components, 3-2 maintaining NEO Iris, 10-1 maintenance, 10-6 as needed tasks, 10-61 adjusting manifold grub screws, 10-84 changing camera lamp, 10-117 checking manifold probes, 10-81 cleaning Liquid Overflow Detection Mat, 10-98 cleaning reader mirror and light diffuser, cleaning the Liquid Level Sensor and Trough, 10-88 emptying system liquid container, 10-63 flushing system liquid, 10-62 performing pipettor verification test, 10-65 performing positions check, 10-67 performing syringe change, 10-68 removing and reinserting manifold, 10-73 removing and replacing probe, 10-105 removing and replacing syringe, 10-99 replacing Y-pusher, 10-123 using the washer teach tool, 10-79 daily tasks, 10-6 checking pipettor reference, 10-17 cleaning instrument, 10-15 completing ABO Anti-A,B reagent quality control, 10-21

completing ABO reagent quality control, 10-20 completing kell phenotyping quality control, 10-22 completing QC3\_Cell quality control, 10-22 deleting historic plate entries, 10-23 empyting common waste container, 10-12 filling system liquid, 10-8 PC and NEO Iris Module Shutdown, 10-7 performing pipettor self check, 10-24 dialog, 3-36 monthly tasks, 10-46 clearing test data, 10-47 decontaminate tubings, 10-49 deleting historic plate entries, 10-48 RVP calibration date, 10-48 records, B-1 weekly tasks, 10-28 checking reader performance, 10-29 checking residual volume, 10-69 cleaning common waste container, 10-43 cleaning pipettor wash towers, 10-40 creating archives, 10-32 inspecting syringes, 10-45 verifying washer, 10-29 maintenance and verification action status accessing, 10-2 maintenance tasks as needed, 10-61 adjusting manifold grub screws, 10-84 changing camera lamp, 10-117 checking manifold probes, 10-81 cleaning Liquid Overflow Detection Mat, 10cleaning reader mirror and light diffuser, cleaning the Liquid Level Sensor and Trough, 10-88 emptying system liquid container, 10-63

flushing system liquid, 10-62 performing pipettor verification test, 10-65 performing positions check, 10-67 performing syringe change, 10-68 removing and reinserting manifold, 10-73 removing and replacing probe, 10-105 removing and replacing syringe, 10-99 replacing Y-pusher, 10-123 using the washer teach tool, 10-79 daily checking pipettor reference, 10-17 cleaning instrument, 10-15 completing ABO Anti-A,B reagent quality control, 10-21 completing ABO reagent quality control, 10-20 completing kell phenotyping quality control, 10-22 completing QC3\_Cell quality control, 10-22 deleting historic plate entries, 10-23 empyting common waste container, 10-12 filling system liquid, 10-8 PC and NEO Iris Module Shutdown, 10-7 performing pipettor self check, 10-24 monthly, 10-46 clearing test data, 10-47 decontaminate tubings, 10-49 deleting historic plate entries, 10-48 RVP calibration date, 10-48 weekly, 10-28 checking reader performance, 10-29 checking residual volume, 10-69 cleaning common waste container, 10-43 cleaning pipettor wash towers, 10-40 creating archives, 10-32 inspecting syringes, 10-45 verifying washer, 10-29 manifold

reinserting, 10-73, 10-78 removing, 10-73 manufacturer information, 1-16 markings, 1-16 microplates barcodes, 1-8 internal barcode scanner, 2-16 technical data, C-5 micro-well strips barcodes, 11

# N

navigation system software, 3-1 Navigation Options, 3-6 Keyboard and Trackball, 3-6 Touch Screen, 3-6 NEO hardware, 2-1 illustration, 2-2 **NEO** Iris extended shutdown, 9-9 instrument start-up, 5-1 intended use, 1-2 logging in, 5-3 logging out, 9-2 maintaining, 10-1 maintenance records, B-1 manufacturer information, 1-16 preparing for first use, A-1 verify all parts are present, A-2 principles of operation, 1-2 repackaging before shipment, A-16 reports, 8-1 safety

user, A-3, A-7 security, 4-1 shutdown, 9-1 shutting down, 9-3 starting up, 5-2 technical support, 1-16 test results, 7-1 troubleshooting, 11-1 using, 6-1 non-crossmatch worklists, II-31, II-32

### 0

on board time plates, 6-32 on/off switch accessibility, A-5

### Ρ

passwords assigning, 4-2 changing, 4-16 PC connections, A-10 descriptions, A-10 personal computer description, 2-3 safety features, 2-4 technical data, C-3 personnel qualified, A-5 pipetting system de-bubbler, 2-23 description, 2-20 pipettors, 2-21

probe wash stations, 2-24 safety features, 2-25 technical data, C-7 pipettor errors recovering from, 11-29 self check failures, 11-14 self-check, 10-24 troubleshooting, 11-27 pipettor reference checking, 10-17 pipettor verification test performing, 10-65 pipettor wash towers cleaning, 10-40 pipettors, 2-21 technical data, C-7 piptest. See pipettor verification test plate based reports accessing, 8-25 body, 8-10 overview, 8-2 plate carrier description, 2-17 improperly loaded, 11-36 improperly taken from centrifuge, 11-38 plate list Assay Protocol tab, 3-30 buttons, 3-24 cancel, 3-27 delete, 3-27 description, 3-23 detail view, 3-29 Flags tab, 3-33

General Information tab, 3-30 information area, 3-23 Plate Events tab, 3-31 plate status list, 3-24 processing steps, 3-26 Raw Results tab, 3-32 Sample ID tab, 3-32 plate loading dialog, 3-64 Assigning Expiration Dates to Plates Expiry Date, 3-69 Automatic Barcode Plate Identification, 3-66 Deselection of Absent Strips Strip Selection, 3-68 Loading Tower Diagram, 3-65 Manual Plate Identification Assay Selection, 3-67 Scan Plate, 3-67 tabs, 3-65 Viewing Processing Steps, 3-72 error locking door, 11-63 error unlocking door, 11-62 technical data, C-4 troubleshooting, 11-62 plate transport troubleshooting, 11-19 plate view, 7-4 plates loading, 6-27 on board time, 6-32 reloading, 6-38 poscheck. See positions check positions check performing, 10-67

post-installation check completing, A-14 power requirements technical data, C-2 principles of operation, 1-2 probe removing, 10-105 replacing, 10-105 probe wash stations, 2-24 process completing sample loading, 6-16 processing starting, 6-34 pump and liquid unit technical data, C-4 pump racks technical data, C-7

# Q

QC reports body, 8-20 QC3\_Cell quality control (QC) completing, 10-22 qualified personnel, A-5

### R

racks, 1-11 parts, 1-11 pump technical data, C-7 types, 1-12 reader checking performance, 10-29 dialog, 3-63 reader mirror and light diffuser cleaning, reagent component grid, I-48 reagent loading bay technical data, C-6 reagents, 1-10 barcodes, 1-8, 10 description, 1-10 loading, 6-21 replenishing, 6-39 reports body, 8-23 reinserting LLS, 10-94 manifold, 10-73, 10-78 reloading plates, 6-38 removing liquid waste, 6-39 LLS, 10-90 manifold, 10-73 probe, 10-105 syringe, 10-99 repackaging NEO Iris, A-16 replacing probe, 10-105 syringe, 10-99 Y-pusher, 10-123 replenishing reagents and controls, 6-39 system liquid, 6-39 reports body, 8-7 current report

body, 8-15 current result overview, 8-2 example, 8-4 footer, 8-8 header, 8-5 NEO Iris, 8-1 overview, 8-2 plate based accessing, 8-25 body, 8-10 overview, 8-2 printing, 8-29 reagents, 8-35 QC body, 8-20 reagents body, 8-23 sample based accessing, 8-27 body, 8-13 overview, 8-2 requesting assays additional, 6-36 requests downloading from LIS, 6-14 requirements sample, 6-3, I-5 rescanning unread barcodes, 6-7 residual volume checking, 10-69 results test interpretation, I-82

viewing, 7-2 RVP verifying, 10-48

### 5

safety user, A-3, A-7 safety features 14-lane and 5-lane bays, 2-19 cabinet, 2-6 camera reader, 2-36 centrifuge, 2-34 hood, 2-10 incubator, 2-27 loading tower, 2-13 personal computer, 2-4 pipetting system, 2-25 transport system, 2-16 washer, 2-30 Safety Precautions, 1-14 sample based reports accessing, 8-27 body, 8-13 overview, 8-2 sample loading process completing, 6-16 STAT Test Processing, 6-16 sample results additional viewing, 7-15 sample view, 7-3 sample/reagenet barcode scanner technical data, C-6 samples

barcodes, 1-9 loading, 6-3, 1-5 additional, 6-36 barcode settings, 6-5 procedure, 6-3 requirements, 6-3, I-5 scanner internal microplate barcode, 2-16 technical data, C-5 sample/reagent barcode technical data, C-6 schedule continual access, 1-14 expected results, 1-14 screen layout, 3-2 machine monitor, 3-3 main menu bar, 3-3 main screen, 3-2 screen layout common buttons, 3-5 confirmation dialogs, 3-5 status bar, 3-3 security NEO Iris, 4-1 selecting assays manually, 6-9 self-check pipettor, 10-24 service, A-5 shutdown after operation, 9-3 dialog, 3-58 extended, 9-9

logging out, 9-2 NEO Iris module, 10-7 PC, 10-7 PC and NEO Iris Module, 10-7 shutting down NEO Iris, 9-3 software configuration, A-13 failures, 11-10 installing, A-12 technical data, C-3 solid waste disposal, 2-38 start run assistant completing the Sample Loading Process, 6-16 description, 3-13 downloading requests from LIS, 6-14 Loading Plates, 6-27 loading Reagents and Controls, 6-21 loading samples, 6-3, I-5 resource overview window, 3-14 Starting Processing, 6-34 steps, 6-2 using, 6-2 starting processing, 6-34 starting up NEO Iris, 5-2 STAT Test Processing, 6-16 Status Bar, 3-81 buttons, 3-81 description, 3-3 Log List, 3-84 Buttons, 3-87 Details Dialog, 3-88 Export Selection Dialog, 3-87

Information Area, 3-85 Log List area, 3-85 Messages, 3-86 Navigation, 3-86 Viewing Messages, 3-84 strips micro-well barcodes, 11 support technical, 1-16 symbols test results, 7-7 syringe inspecting, 10-45 removing, 10-99 replacing, 10-99 syringe change performing, 10-68 syringeex. See syringe change system liquid filling, 10-8 flushing, 10-62 replenishing, 6-39 system liquid container emptying, 10-63 T technical data, C-1 14-lane and 5-lane bays, C-6

cabinet, C-4 camera reader, C-11 centrifuge, C-10 connections, C-2 dimensions, C-2 general, C-2

incubaor, C-9 incubation time, C-9 microplate internal barcode scanner, C-5 personal computer, C-3 pipetting system, C-7 pipettors, C-7 plate loading tower, C-4 power requirements, C-2 pump and liquid unit, C-4 pump rack, C-7 reagent loading bay, C-6 sample/reagent barcode scanner, C-6 software, C-3 transport system, C-4 washer, C-10 weight, C-2 technical support, 1-16 test results interpretation, I-82 test data clearing, 10-47 test details exiting, 7-16 viewing, 7-10 test results approving, 7-17 dialog, 3-34 exiting, 7-19 exporting, 7-18 icons/symbols, 7-7 NEO Iris, 7-1 tool tips, 7-9 view test details, 7-10 viewing, 7-2 event log tab, 7-13
plate view, 7-4 plate views tab, 7-14 reagents tab, 7-12 results tab, 7-12 sample view, 7-3 tool tips, 7-9 **Touch Screen** description, 3-6 Trackball description, 3-6 transport system description, 2-14 internal microplate barcode scanner, 2-16 safety features, 2-16 technical data, C-4 troubleshooting 14-lane and 5-lane bays, 11-60 camera reader, 11-58 centrifuge, 11-32 clot detection, 11-18 error codes, 11-4 incubator, 11-44 NEO Iris, 11-1 pipettor, 11-27 pipettor self check failures, 11-14 plate tower, 11-62 plate transport, 11-19 process steps, 11-2 washer, 11-50 tubings decontaminating, 10-49

## U

unspecified use, A-4 UPS

using, A-11 use unspecified, A-4 user access rights assigning, 4-2 user safety, A-7 users adding, 4-3 changing passwords, 4-16 default access rights, 4-14 deleting, 4-12 editing, 4-10 using NEO Iris, 6-1 start run assistant, 6-2 uninterruptible power supply (UPS), A-11 washer teach tool, 10-79 Utilities dialog, 3-37 Archive tab, 3-38 Event Log tab, 3-38 Printing Statistics, 3-44 Reports tab, 3-45 Statistics tab, 3-39

## V

validating installation, A-15 verifying RVP, 10-48 washer, 10-29 viewing additional sample results, 7-15 archives, 7-20 test results, 7-2 event log tab, 7-13 plate view, 7-4 plate views tab, 7-14 reagents tab, 7-12 results tab, 7-12 sample view, 7-3 viewing test results viewing test details, 7-10 voltage internal, A-4

## W

warnings, 12-16 warranty voiding, A-6 Wash Buffers dialog, 3-62 wash stations probe, 2-24 washer description, 2-28 liquid handling problems, 11-50 mechanical problems, 11-56

module, 2-28 safety features, 2-30 technical data, C-10 troubleshooting, 11-50 using the washer teach tool, 10-79 verifying, 10-29 weight technical data, C-2 work list editor add items, 3-20 add items fields, 3-20 buttons, 3-18 description, 3-17 worklists files crossmatch, II-32 non-crossmatch, II-31, II-32

## Y

Y-pusher replacing, 10-123