## MALB Control

### Package size

1×1 level×1 mL

### Intended use

MALB Control is used in Microalbumin quality control by monitoring precision of Mindray BS measurement system and test ability of laboratory.

### Summary

MALB Control contains Microalbumin (MALB) analyte. Mindray BS measurement system is composed of Mindray BS-series chemistry analyzers, Mindray reagent kits, calibrators and control materials.

### Components

MALB Control is a liquid control based on human albumin. The concentrations of the control components are lot-specific .

### Warnings and Precautions

1. For in vitro diagnostic use only.
2. Professional or trained person is required. Please take the necessary precautions for handling all laboratory reagents.
3. The control values are lot-specific with the matched models listed in the value sheet.
4. The control should be run daily in parallel with the patient samples and after calibration or reagent lot change.
5. Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.
6. The control was tested with Elisa methods and shown to be negative for antibodies to HIV and HCV and nonreactive for HBsAg. However, as no testing method can rule out the potential risk of infection with absolute certainty, so this material should be handled as a patient specimen<sup>1</sup>.
7. Material safety data sheet is available on request.
8. Disposal of all waste material should be in accordance with local guidelines.

### Preparation

1. The control is ready for use. Mix carefully before use. Avoid the formation of foam.
2. Dispense the required volume into a sample cup, and then run QC test you have requested or have it analyzed in the same way as patient samples.

### Storage and Stability

The control is stable up to expiration date indicated on the label when stored in unopened vial at 2~8°C and protected from light. Once opened, it is stable for 14 days when capped tightly and avoiding microbial contamination at 2~8°C.

**Materials required but not provided**

- 1.Mindray BS-series chemistry analyzers, calibrators and reagent kits.
- 2.General laboratory equipment.

**Assay Values**

The control assay values (the target value and range) determined by Mindray standard transfer procedure and routine method are listed in the value sheet. The target value was obtained from Mindray measurement system and the range was calculated as the target values  $\pm 3$  standard deviations. The control values are specific for different lot numbers and different model of chemistry analyzers, so please check the lot number and model of analyzer before use.

**Quality control**

The control result should be within the defined range as shown in the control value sheet. If the control result falls outside the range, the measurement system should be checked. For instance, position of reagent or sample placed on the analyzer; expiration date or storage condition of the calibrator, reagent and control; the parameter setting in software; performance of the analyzer; or the calibration process. Each laboratory should establish its own internal quality control scheme and procedure for corrective action.

**References**

- 1.Biosafety in Microbiological and Biomedical Laboratory;U.S. Department of Health and Human Services; US Government Printing Office;Washington: 2007.

**Graphical symbols**

In Vitro Diagnostic  
medical device

Batch Code

European  
Conformity

Authorized representative in  
the European Community

Use By

Consult  
Instructions for use

Biological  
Risks

Temperature  
Limit

Manufacturer

Control

Catalogue  
number

© 2014 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved

**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

**E-mail Address:** service@mindray.com

**Website:** www.mindray.com

**Tel:** +86-755-81888998

**Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp.  
GmbH(Europe)

**Address:** Eiffestraße 80, Hamburg 20537, Germany

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**Referências**

1. Biosafety in Microbiological and Biomedical Laboratories; U.S. Department of Health and Human Services; US Government Printing Office; Washington: 2007.

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Batch Code

European  
ConformityAuthorized representative in  
the European Community

Use By

Consult  
Instructions for useBiological  
RisksTemperature  
Limit

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**Fabricante:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Endereço:** Mindray building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

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**Üretici:**Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
**AT Temsilcisi:**Shanghai International Holding Corp.GmbH(Europe)