


KNGMED Medical Co. Ltd.	Sodalime TDS Product #: PrSL05-20	Form #: 132 Issued Date: 01.06.2010 Page: 1/1
-------------------------	--------------------------------------	---



TECHNICAL DATA SHEET

SODALIME

Manufacturer Name & Address:	KNGMED Medical Co. Ltd. 10.001 St. No: 28/A-10 Cinar San. Sit. Ulukent, Izmir, TURKIYE 35630	Ph: +902328334262 Fax: +902328334263 Web: www.kngmed.com										
1- Chemical Product Information:	Product Name : Sodalime, (KNGSORB - CO2 absorber – CO2 absorbent) Chemical Name : Sodalime (Medical Carbon dioxide absorbent)											
2- Composition & Ingredients:	Sodium Hydroxide (CAS# 1310-73-2) 2-4% Calcium Hydroxide (CAS# 1305-62-0) >75% Humidity 12-18%											
3- Intended Use	<i>In anesthesia circle systems and respiratory therapy equipment for the purpose of removing exhaled carbon dioxide.</i>											
4- CO₂ Capacity	26% (120 lt) According to USP Methodology											
5- Hardness	75,8 (Shore A)											
6- Particle Distribution	<table><tr><th>Particle Size</th><th>Rate (%)</th></tr><tr><td>> 4 mm</td><td>1</td></tr><tr><td>4-2 mm</td><td>85</td></tr><tr><td>2-0,60mm</td><td>13</td></tr><tr><td>600 µm</td><td>1</td></tr></table>		Particle Size	Rate (%)	> 4 mm	1	4-2 mm	85	2-0,60mm	13	600 µm	1
Particle Size	Rate (%)											
> 4 mm	1											
4-2 mm	85											
2-0,60mm	13											
600 µm	1											
7- Reaction (Indication)	White to Violet											
Please Contact:	KNGMED Medical Co. Ltd. 10.001 St. No: 28/A-10 Cinar San. Sit. Ulukent, Izmir, TURKIYE 35630 Ph: +902328334262 Fax: +902328334263 E-mail: info@kngmed.com & kng.info@yahoo.com											
The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall the company be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential, or exemplary damages howsoever arising, even if the company has been advised of the possibility of such damages.												
Created: On 01.06.2010 by HVK												

Company Profile:	<p>Kngmed Medical was established in 2008 manufactures of whole products in product range. Headquarter is in Izmir, Turkey. It has lots of dealers both Turkey and Abroad.</p> <p>It has 14 employees and capital 2 billion USD.</p>
Executive Summary Product	<p>Sodalime; is a mixture of chemicals, used in granular form in closed breathing environments, absorbs carbon dioxide (CO2) in anesthesia closed circle systems and respiratory therapy equipment.</p> <p>Sodalime; in the close circuits of anesthesia systems, used to absorb CO2 from breathing gases to prevent CO2retention. During the administration of general anaesthesia, the gases expired by a patient, which contain carbon dioxide, are passed through an anaesthetic machine breathing circuit filled with soda lime granules.</p> <p>Production of Sodalime is done newest technology Pelletizer (semi- automatic) and manuel processes. Newly developed automatic systems and technologies are followed and will be applied when the sale quantities and production scale increases.</p> <p>DEVICE DESCRIPTION</p> <p>Product Name: CARBON DIOXIDE ABSORBER (SODALIME)</p> <p>Brand Name: KNGMED</p> <p>Product Models- Photos –</p> 
10.001	<p>SOK NO:28/A -10 ÇINAR SANAYİ SİT. ULUKENT, İZMİR</p> <p>Tel: + 90 232 833 42 62 & Faks: + 90 232 833 42 63</p> <p>www.kngmed.com & kng.info@yahoo.com</p>

<u>Ürünün Adı:</u> <u>Product Name:</u>	<u>Tür/ Boyut:</u> <u>Type/ Size:</u>	<u>Ref No:</u> <u>Ref No:</u>
Sodalime	1 Kg	7600064-1
	5 Kg	7600064-5
	20 Kg	7600064-20

Product MDD Classification and Rule

The product class is determined as **Class II A** according to the **Rule 3 of 93/42/EEC -Annex V** Classification Rules, which is

Class II A Rule 3 Annex V

Medical Device Directive Annex V

The devices are **Class IIa** according to council directive 93/42/EEC, **Annex V, Rule 3** (“This Non-invasive device, modifies chemical composition, which intended for filtrates CO2 breathing circuits of anesthesia devices” classification route is considered).

Biocompatibility Evaluation

N/A Product has no direct or indirect contact with the patient.

Conformity Assessment Method

Medical Devices Directive 93/42/EEC Annex IX .

Product GMDN Code and Description

GMDN Code : **36051**

Category : 02 - Anesthesia and Breathing Devices, 10 – Disposable Devices

Definition : An absorbent material (e.g., granules of treated soda lime) that is placed into an anesthesia system, to remove carbon dioxide (CO2) from the exhaled gases in a patient breathing circuit by chemical reaction. This is a single-use device.

Product Performance Criteria

The performance criteria of our product are determined according to market and user requirements. The requirement controls are based on the 93/42/EEC Medical Devices Directive Annex I Essential Requirements, and Applicable Standards. According to this, determined performance criteria is below:

	<table><tr><th>Criteria</th><th>Claim</th></tr><tr><td>Capacity of CO2 Absorbation</td><td>23% - 28%</td></tr><tr><td>Color Change</td><td>White to Violet</td></tr><tr><td>Grannula Sizes</td><td>4mm R:2mm</td></tr></table>	Criteria	Claim	Capacity of CO2 Absorbation	23% - 28%	Color Change	White to Violet	Grannula Sizes	4mm R:2mm
	Criteria	Claim							
	Capacity of CO2 Absorbation	23% - 28%							
	Color Change	White to Violet							
	Grannula Sizes	4mm R:2mm							
Critical Materials									
Raw Materials / Components									
<table><tr><th>No</th><th>Product Types</th><th>Material</th></tr><tr><td>1</td><td>Chemical Content</td><td>CaOH, NaOH, Ethil Violet</td></tr></table>		No	Product Types	Material	1	Chemical Content	CaOH, NaOH, Ethil Violet		
No	Product Types	Material							
1	Chemical Content	CaOH, NaOH, Ethil Violet							
*: Indicates critical components and raw materials for the product									
Product Storage Information									
Non Sterile products are storage under 0 – 50 °C ± 10 temperature and low humidity (30-70%Rh) conditions before being shipped to the customer. Storage conditions are determined in RP 4.3.2.2 Stability Validation									

Marketing History:	Annual Sale Figures
	2010 Report of Sodalime
	Product Code: Product Name: Sale Amount:
	7600064 Sodalime 6.515
	2011 Report of Sodalime
	Product Code: Product Name: Sale Amount:
	7600064 Sodalime 18.892
	2012 Report of Sodalime
	Product Code: Product Name: Sale Amount:
	7600064 Sodalime 35.986
2013 Report of Sodalime	
Product Code: Product Name: Sale Amount:	
7600064 Sodalime 57.496	
2014 Report of Sodalime	

	Product Code: Product Name: Sale Amount: 7600064 Sodalime 93.984 2015 Report of Sodalime Product Code: Product Name: Sale Amount: 7600064 Sodalime 79.205 2016 Report of Sodalime Product Code: Product Name: Sale Amount: 7600064 Sodalime 76.612 2017 Report of Sodalime Product Code: Product Name: Sale Amount: 7600064 Sodalime 87.541
Work Mechanism:	<p>The product is connected to the anesthesia and ventilation machines. The input gas is flow through the tubes and reach to the patient. After the respiration occurs the gas outlet (expiration) reaches to the machine through the tubes. The gas mixture coming from the patient is directly given into the CO₂ absorber to remove the excess CO₂ content in the expiration gas. The chemical content of the product reacts with the CO₂ and neutralizes its specifications by absorbing the compound.</p> <p>The reaction:</p> <p>1) $H_2O + CO_2 \rightleftharpoons H_2CO_3$ high pH</p> <p>2) $H_2CO_3 + 2 NaOH \rightleftharpoons Na_2CO_3 + 2H_2O$ high pH</p> <p>3) $Na_2CO_3 + Ca(OH)_2 \rightleftharpoons CaCO_3 + 2 NaOH$ high pH also</p> <p>4) $H_2CO_3 + Ca(OH)_2 \rightleftharpoons CaCO_3 + 2H_2O$ high pH</p>
Intended Use :	<p>Sodalime is used for the absorption of carbondioxide gas sourcing from the expiration of the patient who is connected to an anesthesia, ventilation or breathing machine. It doesn't have direct contact with the patient.</p> <p>It is used to provide absorption of CO₂ gas outlet from the patient throughout the anesthesia or ventilation machine.</p>




sodalime

Medical Soda Lime Carbon Dioxide Absorbent

sodalime

(CO₂ absorbent)



Soda lime is essential for CO₂ absorption in inhalation anesthesia machines with rebreathing systems.

With the irregular shape, minimises the amount of dust produced during transit allowing the product to reach you in perfect condition for use.

KNGSORB CO₂ Absorbents (Sodalime) are safely and efficiently removing CO₂ from anesthesia systems, respiratory systems, and hyperbaric chambers.

Characteristic:

- ▶ Particle size: 4x2 mm
- ▶ Hardness: 97%
- ▶ Moisture: 12%-18%
- ▶ CO Absorption: min 130 L/kg
- ▶ CO₂ Activity: min 26%
- ▶ Indicator: White to Violet
- ▶ No KOH

Packages:

- ▶ 5 Kg
- ▶ 20 Kg plastic can

sodalime

(CO₂ absorbent)

KNGSORB is a medical soda lime manufactured from mixture of Calcium and Sodium Hydroxides. It is supplied in the form of hard, porous regular rounded pellets which have been specifically processed to absorption CO₂ capacity. One kilogram of KNGSORB will absorb a minimum of 130 litres of CO₂.



KNGSORB does not contain potassium Hydroxide or Barium Hydroxide and only low level of Sodium Hydroxide. High CO₂ absorption capacity pellet form prevents channel formation, uniform utilization reducing your costs through longer life time. Little dust generation creating a safer environment for your personnel and protects your equipment. Additional safety for patients. Reduction of undesirable reaction products when accidentally dried out.

Raw material

- ▶ NaOH less than 4%
- ▶ KOH free formula
- ▶ 14-19% H₂O

Design

- ▶ Optimum air flow and absorption 4x2 mm granule form

Features

- ▶ 97% hardness
- ▶ Color indicator shows usable life changes from white to violet

Usage

- ▶ Suitable for all close circuit anesthesia devices

Package

- ▶ 5kg and 20kg plastic can
- ▶ Light proof, opaque package
- ▶ Leakage proof, sealed lid



Storage information

Store the absorbent in original closed container only. Use within one month of opening the container. Keep container away from direct heat and sunlight to prevent drying out or over-heating of the absorbent. Do not freeze; if frozen, absorbent must be completely thawed before use. Do not subject to excessive loads. When stored correctly, unopened containers will maintain their absorption capability for at least five years.



Safety information

Risk of serious damage to eyes. Irritating to eye, respiratory system and skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable gloves and eye/face protection.

1. Sodalime 5kg - item no: #7600064-5
2. Sodalime 20kg (item no: #7600064-20)