Medtronic

Instructions for use

React™ Catheter

fr	Cathéter
de	Katheter
it	Catetere
es	Catéter
sv	Kateter
nl	Katheter
pt	Cateter
fi	Katetri
da	Kateter
el	Καθετήρας
cs	Katétr
hu	Katéter
ru	Катетер
pl	Cewnik
tr	Kateter
no	Kateter
sk	Katéter
ro	Cateter
bg	Катетър
et	Kateeter
hr	Kateter
lt	Kateteris
lv	Katetrs
sl	Kateter
sr	Kateter
pt-br	Cateter
ko	카테터
vi	Ống thông
kk	Катетер
mk	Катетер
id	Kateter
ar	قسطرة

TABLE OF CONTENTS

React™ Catheter

English	3
Français	4
Deutsch	5
Italiano	6
Español	7
Svenska	8
Nederlands	9
Português	10
Suomi	11
Dansk	12
Ελληνικά	13
Česky	15
Magyar	16
Русский	17
Polski	18
Türkçe	20
Norsk	21
Slovenčina	22
Română	23
Български	24
Eesti	25
hrvatski	26
Lietuvių	27
Latviski	28
Slovenščina	29
Srpski	30
Português (Brasil)	31
한국어	32
Tiếng Việt	34
Қазақша	35
Македонски	36
Bahasa Indonesia	38
العربية	39
Symbol Glossary	40

Instructions for Use

React™ Catheter

CAUTION

- Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician.
- This device should be used only by physicians trained in percutaneous, intravascular techniques and procedures.

DESCRIPTION

The React™ Catheter is a single lumen, flexible, variable stiffness composite catheter. The catheter shaft has a hydrophilic coating that spans the distal 40 cm to reduce friction during use. The catheter shaft is visible under fluoroscopy. The catheter dimensions are included on the individual device label. The catheter is introduced through a guide catheter or sheath and into the intracranial vasculature and guided over a neurovascular guidewire and/or micro catheter. The proximal end of the catheter has a luer fitting to allow attachment of accessories and infusion of liquids through the system.


Figure 1. React™ Catheter

- Effective length
- Distal tip outer diameter
- Distal inner diameter
- Proximal shaft outer diameter
- Proximal inner diameter

Table 1. Dimensions of React™ Catheter

CFN	Effective length (1)	Distal tip outer diameter (2)	Distal inner diameter (3)	Proximal shaft outer diameter (4)	Proximal inner diameter (5)	Guidewire Maximum Diameter
REACT-68	132 cm	6.3 F 0.083" (2.11 mm)	0.068" (1.73 mm)	6.3 F 0.083" (2.11 mm)	0.068" (1.73 mm)	0.038" (0.97 mm)
REACT-71	132 cm	6.5 F 0.0855" (2.17 mm)	0.071" (1.80 mm)	6.5 F 0.0855" (2.17 mm)	0.071" (1.80 mm)	0.038" (0.97 mm)

DEVICE COMPATIBILITY

The maximum diameter recommended for a guidewire is 0.038" (0.97 mm).

When using the React™ 71 Catheter, the minimum inner diameter recommended for a sheath is 0.087" (2.21 mm).

Refer to product label and Table 1 for device dimensions. Refer to labeling provided with other medical technologies to determine compatibility.

INTENDED PURPOSE / INDICATIONS FOR USE

The React™ 68 Catheter and React™ 71 Catheter are indicated for the introduction of interventional/diagnostic devices into the peripheral and neurovasculature.

The React™ 68 Catheter and React™ 71 Catheter are also indicated for the removal/aspiration of emboli and thrombi from selected blood vessels in the arterial system, including the neurovasculature.

CONTRAINDICATIONS

There are no known contraindications.

PREPARATIONS FOR USE

- Gently remove the catheter and packaging hoop from the pouch by grasping the packaging hoop and slowly pulling the catheter and hoop out of the pouch.
- Flush package hoop with saline before removing catheter.
- Remove the catheter from the packaging hoop.
- Inspect the catheter for kinks or other damage. If any damage is observed, replace with new device.
- Connect a rotating hemostasis valve to the hub of the catheter and flush the inner lumen with heparinized saline.

DIRECTIONS FOR USE

- Introduce the catheter into the vasculature through the sheath, and over a guide wire using a percutaneous entry technique of choice. A split sheath is provided in the packaging to provide support and facilitate the introduction of the tip of the catheter into the sheath.
NOTE: The maximum diameter recommended for a guide wire is 0.038" (0.97 mm). When using the React™ 71 Catheter, the minimum inner diameter recommended for a sheath is 0.087" (2.21 mm).
- Under fluoroscopic guidance, advance the system over the guide wire until the desired position is attained.
- Remove the guide wire.
- If introducing interventional devices, refer to the instructions for use provided by the manufacturer.
- If performing aspiration, apply aspiration to the catheter.
- If necessary, additional attempts at flow restoration may be used.
- Obtain a post-treatment angiogram in accordance with the physician's preferred technique.

POTENTIAL COMPLICATIONS

Potential complications of the device and the endovascular procedure include or are synonymous with, but may not be limited to the following:

- Access site complications such as hematoma, pain, nerve damage
- Adverse reactions to antiplatelet/anticoagulation agents, contrast media, or anesthesia such as hypersensitivity, hypotension, hemorrhage
- Cardiac complications such as arrhythmia
- Complications of radiation exposure such as alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia
- Device complications such as friction, foreign body in patient, device separation, device fracture
- Hematologic complications such as coagulopathy, hemolysis, intracranial hemorrhage
- Neurological deficit/dysfunctions such as cerebral infarction, transient ischemic attack, cognitive changes, hydrocephalus, edema, foreign body reaction, organ failure, toxicity, shock
- Systemic complications such as infection, fever, inflammation
- Therapeutic response decreased
- Vascular complications such as dissection, perforation, rupture, ischemia, embolism (including distal embolization to previously uninvolved territory), thromboembolism, thrombosis, fistula, necrosis, stenosis, vasoconstriction, occlusion
- Vision symptoms such as impaired vision and vision loss
- Death

* Consult instructions for use for other therapy devices and medications for additional potential complication information. If a serious incident related to the device occurs, contact your Medtronic representative and the competent authority in your respective country/region.

WARNINGS

- The catheter should only be used by physicians who have received appropriate training in interventional techniques.
- Do not reuse. The device is intended for single use only. Discard the catheter after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- The catheter has not been tested for use with automated high-pressure contrast injection equipment, do not use this equipment with the device because it may damage the device.
- Visually inspect all sterile barrier systems, that are labeled as sterile, immediately prior to use. Do not use the device if breaches in sterile barrier system integrity are evident.
- For additional Materials of Concerns information such as REACH, CA Prop 65 or other product stewardship programs, go to www.medtronic.com/productstewardship

PRECAUTIONS

- Do not use open, kinked, or damaged devices.
- Do not autoclave.
- See the product label for the device "Use-by date." Do not use the device beyond its labeled "Use-by date."
- Use the catheter in conjunction with fluoroscopic visualization.
- Inspect the catheter before use to verify that its size and condition are suitable for the specific procedure.
- Do not advance or withdraw the catheter against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Moving or torquing the device against resistance may result in damage to the vessel or device.
- Maintain a constant infusion of appropriate flush solution.
- If flow through the catheter becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.
- Extreme care must be taken to avoid damage to the vasculature through which the catheter passes. The catheter may occlude smaller vessels. Care must be taken to avoid complete blood flow blockage.
- Torquing the catheter may cause damage which could result in kinking and possible separation along the catheter shaft. Should the catheter become severely kinked, withdraw the catheter.
- An appropriate anticoagulation therapy should be applied per institutional guidelines.
- Operators should take all necessary precautions to limit X-radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

HOW SUPPLIED

- This device is supplied STERILE using ethylene oxide. This device is Non-pyrogenic.

STORAGE AND DISPOSAL

- This device should be stored in a dry place, away from sunlight.
- Dispose of device in accordance with hospital, administrative, and/or local government policy.

SAFETY AND CLINICAL PERFORMANCE

- Expected clinical benefit - The React™ 68 Catheter and React™ 71 Catheter are used to access the peripheral and neurovasculature to attain the following clinical benefits:
 - Introduction of interventional/diagnostic devices
 - Removal/aspiration of emboli and thrombi from selected blood vessels
- Performance characteristics of the device - The React™ Catheter is a single lumen, flexible, variable stiffness composite catheter. The React™ Catheter introduces interventional/diagnostic devices into the peripheral and neurovasculature, and removes/aspirates emboli and thrombi from selected blood vessels in the arterial system, including the neurovasculature. A copy of the Summary of Safety and Clinical Performance (SSCP) can be viewed by searching the device name on the EUDAMED website <https://ec.europa.eu/tools/eudamed>. Search for the SSCP using the manufacturer and device name, and any of the following elements, as applicable: device model, reference number, catalog number, or the Basic Unique Device Identification (Basic UDI-DI) number-847536080000165J.

DISCLAIMER WARRANTY

Although this product has been manufactured under carefully controlled conditions, the manufacturer has no control over the conditions under which this product is used. The manufacturer therefore disclaims all warranties, both expressed and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. The manufacturer shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind the manufacturer to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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© 2023 Medtronic.
Medtronic, Medtronic with logo and Medtronic logo are U.S. and internationally registered trademarks of Medtronic.
Micro Therapeutics, Inc.
Medtronic B.V.