DCM001-12  
Ed. 11/2016

# CORTISOL ELISA

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# CORTISOL ELISA

for routine analysis

Direct immunoenzymatic determination of Cortisol in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO001

## INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Cortisol concentration in human serum or plasma.

Cortisol ELISA kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

Cortisol is a steroid hormone released from the adrenal cortex in response to a hormone called ACTH (produced by the pituitary gland), it is involved in the response to stress; it increases blood pressure, blood sugar levels, may cause infertility in women, and suppresses the immune system.

Cortisol acts through specific intracellular receptors and has effects in numerous physiologic systems, including immune function, glucose-counter regulation, vascular tone, substrate utilization and bone metabolism. Cortisol is excreted primarily in urine in an unbound (free) form.

Cortisol is bound, in plasma, from corticosteroid-binding globulin (CBG, transcortin), with high affinity, and from albumin. Only free cortisol is available to most receptors.

The amount of cortisol present in the serum undergoes diurnal variation, with the highest levels present in the early morning, and lower levels in the evening, several hours after the onset of sleep. Highest levels are at about 6-8 a.m. and lowest levels are at about midnight. These normal endogenous functions are the basis for the physiological consequences of chronic stress - prolonged cortisol secretion causes muscle wastage, hyperglycaemia, and suppresses immune / inflammatory responses. The same consequences arise from long-term use of glucocorticoid drugs.

## 2. PRINCIPLE OF THE METHOD

The Cortisol (antigen) in the sample competes with the antigenic Cortisol conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti Cortisol coated on the microplate (solid phase).

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Then, the enzyme HRP in the bound-fraction reacts with the Substrate (H<sub>2</sub>O<sub>2</sub>) and the TMB Substrate and develops a blue color that changes into yellow when the Stop Solution (H<sub>2</sub>SO<sub>4</sub>) is added.

The colour intensity is inversely proportional to the Cortisol concentration of in the sample.

Cortisol concentration in the sample is calculated through a calibration curve.

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REF DCE002/0110-0

2. Control (1 vial, 1 mL)

See Control concentration on the Certificate of Analysis

REF DCE045/0103-0

3. Conjugate (1 vial, 21 mL)

Cortisol conjugated with Horseradish peroxidase (HRP)

REF DCE002/0102-0

4. Coated Microplate (1 breakable microplate)

Microplate coated with anti Cortisol antibodies

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## 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of Cortisol from 10 ng/mL to 500 ng/mL.
- The clinical significance of the Cortisol determination can be invalidated if the patient was treated with corticosteroids or natural or synthetic steroids.

## 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.

- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibration (C<sub>0</sub>...C<sub>4</sub>)

The Calibrators are ready for use and have the following concentration of Cortisol:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>
ng/mL	0	10	50	150	500

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### 6.3. Preparation of the Sample

The determination of Cortisol with this kit can be performed in human plasma as well as in serum.

Store the sample at -20°C if the determination is not performed on the same day of the sample connection. Avoid repetitive freezing and thawing of samples.

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In concentrated wash solution is possible to observe the presence of crystals; in this case mix at room

temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

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Reagent	Calibrator	Sample/ Control	Blank
Calibrator C <sub>0</sub> -C <sub>4</sub>	20 µL		
Sample/ Control		20 µL	
Conjugate	200 µL	200 µL	
Incubate 1 hour at 37°C. Remove the contents from each well. Wash the wells 3 times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> in case you use an automatic washer, it is advised to do 6 washing steps.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22÷28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Cortisol for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

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### 8.2. Calibration Curve

Plot the mean value of absorbance (Em) of the calibrators (C<sub>0</sub>-C<sub>4</sub>) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

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The serum or plasma Cortisol reference values are:

60 - 230 ng/mL between 8.00 – 10.00 A.M.  
 30 – 150 ng/mL at 4.00 P.M.

Patient treated with ACTH: 280 - 600 ng/mL  
 Patient treated with dexamethasone: 0 - 50 ng/mL

## 10. PERFORMANCE AND CHARACTERISTICS

### 10.1. Precision

#### 10.1.1. Intra Assay

Within run variation was determined by replicate (16x) the measurement of three different sera in one assay. The within assay variability is 5.1%.

#### 10.1.2. Inter Assay

Between run variation was determined by replicate (10x) the measurement of three different sera in different lots. The between assay variability is 11.0%.

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### 10.3. Specificity: interfering substances

Interference by Bilirubin, Hemoglobin and Triglycerides has been investigated in Diametra Cortisol ELISA kit:

Substance	Assayed Conc.	Interference
Bilirubina	0.2 mg/mL	No
Emoglobina	2 mg/mL	No
Trigliceridi	6 mg/mL	No

No interference has been observed with the substances under investigation; anyway, according to good laboratory practices, avoid to use highly lipemic or haemolysed samples.

### 10.4. Specificity: plasma and SST tube

Interference in plasma and SST (serum separation tube) samples has been evaluated. Serum obtained from the same patient has been used as reference.

Sample	Interference
SST (serum separation tube)	No
EDTA plasma	No
Litio-eparina plasma	No
Sodio-eparina plasma	No

No interference has been observed.

### 10.5. Accuracy

The recovery of 15 - 35 - 70 ng/mL of Cortisol added to samples gave an average value ( $\pm$ SD) of 101.61%  $\pm$  4.42% with reference to the original concentrations. The dilution test performed on three sera diluted 2 - 4 and 8 times gave an average value ( $\pm$ SD) of 106.18%  $\pm$  4.68%.

### 10.6. Sensitivity

The lowest detectable concentration of cortisol that can be distinguished from the Calibrator zero is 2.42 ng/mL at the 95% confidence limit.

### 10.7. Correlation

Diametra Cortisol ELISA kit has been compared to a commercially available chemiluminescence Cortisol kit. 61 serum samples have been analysed. The linear regression curve has been calculated:  
 $Y = 1.09 \cdot X - 19.23$   
 $r^2 = 0.904$

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**ERROR POSSIBLE CAUSES / SUGGESTIONS****No colorimetric reaction**

- no conjugate pipetted reaction after addition
- contamination of conjugates and/or of substrate
- errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)

**Too low reaction (too low ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too short, incubation temperature too low

**Too high reaction (too high ODs)**

- incorrect conjugate (e.g. not from original kit)
- incubation time too long, incubation temperature too high
- water quality for wash buffer insufficient (low grade of deionization)
- insufficient washing (conjugates not properly removed)

**Unexplainable outliers**

- contamination of pipettes, tips or containers
- insufficient washing (conjugates not properly removed) too high within-run
- reagents and/or strips not pre-warmed to CV% Room Temperature prior to use
- plate washer is not washing correctly (suggestion: clean washer head)
- too high between-run - incubation conditions not constant (time, CV % temperature)
- controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)
- person-related variation

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# TESTOSTERONE ELISA

for routine analysis

Direct immunoenzymatic determination of total Testosterone in human serum or plasma.

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO002

## INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of total Testosterone concentration in human serum or plasma.

Testosterone ELISA kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

Testosterone (17β-OH-4-androstene-3-one) is a steroid hormone family of androgens mainly produced by the Leydig cells located in the testes and, minimally, by the ovaries and the adrenal cortex. It is also present in women who, compared to men, have a greater tendency to convert into estrogen this hormone.

In postpubertal males, testosterone is secreted primarily by the testes with only a small amount derived from peripheral conversion of androstenedione. In humans is deputy to the development of the sexual organs (differentiation of the testis and the whole genital apparatus) and of secondary sexual characteristics, such as beard, hair distribution, the timbre of the voice and muscles. The testosterone, during puberty, is also involved on skeletal development, limiting the elongation of the long bones and avoiding, in this way, a disproportionate growth of the limbs. In adult humans, the levels of testosterone have a very important role as regards the sexuality, the musculoskeletal system, the vitality and good health (mainly understood as protection from metabolic diseases such as hypertension and diabetes mellitus); helps to ensure fertility, as it stimulates the maturation of sperm in the testes. Also influence the quality and quantity of sperm produced, since the seminal work on the streets and on the prostate, deputies to the production of sperm. Daily production of testosterone in men varies from 5 to 7 milligrams, but exceeded 30 years, tends to decrease annually by 1%. In adult women over 50% of serum testosterone is derived from peripheral conversion of androstenedione secreted by the adrenal and ovary, with the remainder from direct secretion of testosterone by these glands. The majority of circulating testosterone is bound by SHBG and a smaller portion is bound by albumin. Only a small percentage (< 1%) exists in circulation as unbound or free testosterone.

Testosterone effects can be classified as virilizing and anabolic effects, although the distinction is somewhat artificial, as many of the effects can be considered both. Anabolic effects include growth of muscle mass and strength, increased bone density and strength, and stimulation of linear growth and bone maturation. Virilizing effects include maturation of the sex organs, and after birth (usually at puberty) a deepening of the voice, growth of the beard and axillary hair (male secondary sex characteristics).

Testosterone levels decline gradually with age in men (andropause). The signs and symptoms are non-specific, and are generally associated with aging such as loss of muscle mass and bone density, decreased physical endurance, decreased memory ability and loss of libido.

In females of all ages, elevated testosterone levels can be associated with a variety of virilizing conditions, including adrenal tumors and polycystic ovarian disease.

## 2. PRINCIPLE

Testosterone in the blood is bound to SHBG (60%) and in lower quantity to other proteins (for example albumin); the unbound Testosterone (< 1% of total Testosterone) is known as "free Testosterone". The chemical formulation of this assay allows to release completely the Testosterone from bound proteins; thus Diametra Testosterone kit allows to measure the concentration of total Testosterone (bound + free) in the sample. For the measurement of free Testosterone only, ELISA Diametra "Free Testosterone" kit is available.

Testosterone (antigen) in the sample competes with the antigenic Testosterone conjugated with horseradish peroxidase (HRP) present in the Conjugate for binding to the antibodies anti-testosterone coated on the microplates (solid phase). After the incubation, the bound/free separation is performed by a simple solid-phase washing.

The enzyme HRP in the bound-fraction reacts with the Substrate (H<sub>2</sub>O<sub>2</sub>) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H<sub>2</sub>SO<sub>4</sub>) is added.

The colour intensity is inversely proportional to the Testosterone concentration in the sample.

Testosterone concentration in the sample is calculated through a calibration curve.

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Controls Concentration is indicated on the Certificate of Analysis

##### 3. Conjugate (1 vial, 12 mL)

Testosterone conjugated with horseradish peroxidase (HRP)

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Anti Testosterone antibody adsorbed on microplate

**REF** **DCE002/0203-0**

##### 5. TMB Substrate (1 vial, 15 mL)

H<sub>2</sub>O<sub>2</sub>-TMB 0.26 g/L (avoid any skin contact)

**REF** **DCE004-0**

##### 6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)

**REF** **DCE005-0**

##### 7. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M pH 7.4

**REF** **DCE054-0**

#### 3.2. Reagents necessary not supplied

Distilled water

#### 3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

#### Notes

Store all reagents at 2-8 °C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, the microplate is stable until expiry date of the kit.

#### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagents should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as

healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.

- Some reagents contain small amounts of Sodium Azide (NaN<sub>3</sub>) or Proclin 300<sup>R</sup> as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of Testosterone from 0.2 ng/mL to 16.0 ng/mL.
- The clinical significance of the determination Testosterone can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.

#### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8 °C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28 °C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems, it is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the

addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>4</sub>)

The Calibrators are ready for use and have the following concentration of Testosterone:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>
ng/mL	0	0.2	1.0	4.0	16.0

Once opened, the Calibrators are stable six months at 2-8°C.

### 6.2. Preparation of the Conjugate

Ready to use. Mix gently for 5 minutes with a rotating mixer.

Once opened, it is stable six months at 2-8°C.

### 6.3. Preparation of the Sample

The determination of Testosterone can be performed in human plasma as well as in serum of patients.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Avoid repetitive freezing and thawing of samples.

The Controls are ready to use.

### 6.4. Preparation of Wash Solution

Dilute the content of each vial of the "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

In concentrated wash solution is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

### 6.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.

- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>4</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Sample/Control		25 µL	
Calibrator C <sub>0</sub> -C <sub>4</sub>	25 µL		
Conjugate	100 µL	100 µL	
Incubate at 37°C for 1 hour. Remove the content from each well. Wash the wells 3 times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> in case you use an automatic washer, it is advised to do 6 washing steps.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22÷28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Testosterone for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 8. RESULTS

### 8.1. Mean absorbance

Calculate the mean of the absorbances ( $E_m$ ) for each point ( $C_0$ - $C_4$ ) of the calibration curve and of each sample.

### 8.2 Calibration curve

Plot the mean value of absorbance ( $E_m$ ) of the Calibrators ( $C_0$ - $C_4$ ) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

### 8.3 Calculation of results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in ng/mL.

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### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate the measurements (16x) of three different sera in one assay. The within assay variability is  $\leq 7.0\%$ .

#### 10.1.2. Inter Assay Variation

Between run variation was determined by replicate the measurements (9x) of three different sera in different lots. The between assay variability is  $\leq 8.3\%$ .

### 10.2. Specificity

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

Testosterone	100%
Dihydrotestosterone	2.03%
Androstenedione	0.01%
Androsterone	0.05%
DHEA-S	0%
Cortisol	0.01%
Cortisone	0%
17b Estradiol	0.16%
Estrone	0%
Prednisone	0.01%

### 10.3. Accuracy

The recovery of 0.4 - 0.8 - 4.0 - 14.0 ng/mL of Testosterone added to the sample gave an average value of 98.9% with reference to the original concentrations.

The dilution test performed on 3 samples diluted up to 4 times gave an average value of 99.7%.

### 10.4. Sensitivity

The lowest detectable concentration of Testosterone that can be distinguished from the Calibrator 0 is 0.10 ng/mL at the 95% confidence limit.

### 10.5. Correlation

The new Diametra Testosterone ELISA kit was compared to the old Diametra Testosterone ELISA kit. Serum samples from 24 females and 28 males were analysed.

The linear regression curve was calculated:

$$Y = 1.13 \cdot X + 0.06$$

$$r^2 = 0.92$$

Diametra Testosterone ELISA kit was compared to another commercially available Testosterone assay. Serum samples of 24 females and 29 males were analysed according in both test systems.

The linear regression curve was calculated:

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Tel. +39-02-2139184

Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14,  
06038 SPELLO (PG) Italy

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REF DKO003

## INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of 17 $\beta$ -Estradiol concentration in human serum or plasma.

Estradiol ELISA kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

Estradiol (17 $\beta$ -Estradiol) is a sex hormone. It represents the major estrogen in humans. Estradiol has not only a critical impact on reproductive and sexual functioning, but also affects other organs including bone structure. During the reproductive years most Estradiol in women is produced by the ovaries, smaller amounts of Estradiol are also produced by the adrenal cortex. In men, the testes produce Estradiol. In plasma Estradiol is largely bound to sex hormone binding globulin (SHBG), also to albumin, only a fraction is free and biologically active. Serum Estradiol measurement in women reflect primarily the activity of the ovaries. During pregnancy estrogen levels, including Estradiol, rise steadily towards term. Estradiol increases due to placental production. In adult premenopausal women, ovarian production of Estradiol is stimulated by luteinizing hormone (LH) and the follicle-stimulating hormone (FSH) during the menstrual cycle.

In adult women, Estradiol levels are measured in the evaluation of fertility and menstrual irregularities, and to monitor ovarian follicular function during induction of ovulation. In the female, Estradiol acts as a growth hormone for tissue of the reproductive organs. The development of secondary sexual characteristics in women is driven by Estradiol. Estradiol is involved also in man fertility. Estradiol regulates the bone maintenance. Post-menopause women experience an accelerated loss of bone mass due to a relative estrogen deficiency.

Estradiol affects the production of multiple proteins including lipoproteins, binding proteins, and proteins responsible for blood clotting. Estrogens have been found to have neuroprotective function. The Estradiol, for his activities, is involved in some types of cancer such as breast cancer and cancer of the uterine lining. In addition there are several benign gynecologic conditions that are dependent on

estrogen such as endometriosis, leiomyomata uteri, and uterine bleeding.

## 2. PRINCIPLE

17 $\beta$ -Estradiol (antigen) in the sample competes with the antigenic 17 $\beta$ -Estradiol conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti 17 $\beta$ -Estradiol coated on the microplate (solid phase).

After incubation, the bound/free separation is performed by a simple solid-phase washing. Then the enzyme HRP in the bound-fraction reacts with the Substrate (H<sub>2</sub>O<sub>2</sub>) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H<sub>2</sub>SO<sub>4</sub>) is added. The colour intensity is inversely proportional to the 17 $\beta$ -Estradiol concentration in the sample.

17 $\beta$ -Estradiol concentration in the sample is calculated through a calibration curve.

## 3. REAGENT, MATERIAL AND INSTRUMENTATION

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17 $\beta$ -Estradiol conjugated with horseradish peroxidase (HRP)

REF DCE002/0302-0

#### 4. Coated Microplate (1 microplate breakable)

Anti 17 $\beta$ -Estradiol antibodies adsorbed on microplate

REF DCE002/0303-0

#### 5. TMB Substrate (1 vial, 15 mL)

H<sub>2</sub>O<sub>2</sub>-TMB (0.26 g/L) (avoid any skin contact)

REF DCE004-0

#### 6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

7. 10X Conc. Wash Solution (1 vial, 50 mL)  
Phosphate buffer 0.2M pH 7.4 **REF** DCE054-0

### 3.2. Reagents necessary not supplied in the kit

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Microplates reader (450 nm, 620-630 nm)

### Notes

*Store all reagents between 2÷8 °C in the dark.*

*Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.*

### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of 17β-Estradiol from 20 pg/mL to 2000 pg/mL.
- The clinical significance of 17β-Estradiol determination can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.

### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8 °C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28 °C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.

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### 6. PROCEDURE

#### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are ready to use and have the following concentration of 17β-Estradiol:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
pg/mL	0	20	120	300	600	2000

The Calibrators are stable until the expiry date printed on the label. Once opened, the standards are stable six months at 2÷8 °C.

#### 6.2. Preparation of the Conjugate

The conjugate is ready to use. Mix gently, for 5 minutes, with a rotating mixer.

Once opened, it is stable six months at 2÷8 °C

#### 6.3. Preparation of Wash Solution

Dilute the content of the vial "Conc. Wash Solution 10X" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2÷8 °C. In concentrated wash solution is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals, for greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL taking care also to transfer crystals, then mix until crystals are completely dissolved.

#### 6.4. Preparation of the sample

The determination of 17 $\beta$ -Estradiol can be performed in human serum or plasma. Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Avoid repetitive freezing and thawing of samples. Before using, mix gently, for 5 minutes, with a rotating mixer. The Control is ready for use.

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Reagent	Calibrator	Sample/ Control	Blank
Calibrator C <sub>0</sub> -C <sub>5</sub>	25 $\mu$ L		
Sample/ Control		25 $\mu$ L	
Conjugate	200 $\mu$ L	200 $\mu$ L	
Incubate 2 h at +37°C. Remove the contents from each well, wash the wells three times with 300 $\mu$ L of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 $\mu$ L	100 $\mu$ L	100 $\mu$ L
Incubate at room temperature (22 $\div$ 28°C) for 30 minutes in the dark.			
Stop Solution	100 $\mu$ L	100 $\mu$ L	100 $\mu$ L
Shake gently the microplate. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

#### 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Estradiol for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

#### 8. RESULTS

##### 8.1. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>) and of each sample

##### 8.2. Calibration curve

Plot the mean value of absorbance (Em) of the calibrators (C<sub>0</sub>-C<sub>5</sub>) against concentration. Draw the best-fit curve through the plotted points (es: Four Parameter Logistic).

##### 8.3. Calculation of results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in pg/mL.

#### 9. REFERENCE VALUES

The serum 17 $\beta$ -Estradiol reference values are:

		pg/mL
<b>WOMEN</b>	Follicular phase	30 – 100
	Ovulatory peak	130 – 350
	Luteinic phase	50 – 180
	Menopause	< 60
<b>CHILDREN</b>		< 40
<b>MEN</b>		< 60

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**10.1. Precision***10.1.1. Intra Assay Variation*

Within run variation was determined by replicate measurements (10x) of three different human sera in one assay. The within assay variability is  $\leq 9\%$ .

*10.1.2. Inter Assay Variation*

Between run variation was determined by replicate measurements of three different human sera in different lots. The between assay variability is  $\leq 10\%$ .

**10.2. Accuracy**

The dilution test conducted with high concentration samples of  $17\beta$ -Estradiol gave an average recovery value ( $\pm$ SD) of  $95.69\% \pm 7.74\%$  with reference to the original concentrations.

The recovery of 120 – 240 – 480 – 960 pg/mL of  $17\beta$ -Estradiol added to samples gave an average value ( $\pm$ SD) of  $101.09\% \pm 5.42\%$  with reference to the original concentrations.

**10.3. Sensitivity**

The lowest detectable concentration of  $17\beta$ -Estradiol that can be distinguished from the calibrator zero is 8.68 pg/mL at the 95% confidence limit.

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**10.5. Correlation with a commercial kit**

Diametra Estradiol ELISA kit was compared to another commercially available  $17\beta$ -Estradiol assay. 16 serum samples were analysed in both test systems.

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# PROGESTERONE ELISA

for routine analysis


Direct immunoenzymatic determination of Progesterone in human serum or plasma

IVD

LOT

See external label

$\Sigma = 96$  tests

REF DKO006

## INTENDED USE

Competitive immunoenzymatic colorimetric method for quantitative determination of Progesterone concentration in human serum or plasma.

Progesterone ELISA kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

Progesterone is a C-21 steroid hormone involved in the female menstrual cycle, pregnancy (supports *gestation*) and embryogenesis of humans and other species. Progesterone is the major naturally occurring human progestagen.

Progesterone is important for aldosterone (mineralocorticoid) synthesis, as 17-hydroxyprogesterone is for cortisol (glucocorticoid).

Progesterone levels are relatively low in children and postmenopausal women. Adult males have levels similar to those in women during the follicular phase of the menstrual cycle.

In women, progesterone levels are relatively low during the preovulatory phase of the menstrual cycle, rise after ovulation, and are elevated during the luteal phase. If pregnancy occurs, progesterone levels are maintained at luteal levels initially. After delivery of the placenta and during lactation, progesterone levels are very low. The fall in progesterone levels following delivery is one of the triggers for milk production.

Progesterone is produced in the adrenal glands, the gonads (specifically after ovulation in the corpus luteum), the brain, and, during pregnancy, in the placenta.

Progesterone converts the endometrium to its secretory stage to prepare the uterus for implantation. If pregnancy does not occur, progesterone levels will decrease, leading, in the human, to menstruation.

Progesterone belongs to the group of neurosteroids that are found in high concentrations in certain areas in the brain and are synthesized there. Neurosteroids affect synaptic functioning, are neuroprotective, and affect myelinization.

Progesterone has multiple effects outside of the reproductive system. Progesterone is thermogenic, it reduces spasm and relaxes smooth muscle. Bronchi are widened and mucus regulated. Progesterone acts as an antiinflammatory agent and regulates the immune response. Progesterone also assists in thyroid function, in bone building by osteoblasts.

Measurement of serum progesterone concentrations have been used in evaluating ovarian function.

## 2. PRINCIPLE

The Progesterone (antigen) in the sample competes with the antigenic Progesterone conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti Progesterone coated on the microplate (solid phase).

After incubation, the bound/free separation is performed by a simple solid-phase washing.

Then, the enzyme HRP in the bound-fraction reacts with the Substrate ( $H_2O_2$ ) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution ( $H_2SO_4$ ) is added.

The colour intensity is inversely proportional to the Progesterone concentration of in the sample.

Progesterone concentration in the sample is calculated through a calibration curve.

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REF DCE045/0603-0

3. Conjugate (1 vial, 22 mL)

Progesterone conjugated with horseradish peroxidase (HRP)

REF DCE002/0602-0

4. Coated Microplate (1 microplate breakable)

Antibody anti Progesterone coated on the microplate

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## 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants.
- This method allows the determination of Progesterone from 0.2 ng/mL to 40.0 ng/mL.
- For higher values, for example in pregnancy, dilute the sample; consider the diluting factor when calculating the final result.
- The clinical significance of the Progesterone determination can be invalidated if the patient was treated with cortisone or natural or synthetic steroids.

## 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.

- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrator (C<sub>0</sub>...C<sub>4</sub>)

Before using leave 5 minutes on a rotary shaker.

The Calibrators are ready to use and have the following concentrations of Progesterone:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>
ng/mL	0	0.2	1.0	8.0	40.0

The Calibrators are stable until the expiry date printed on the label. Once opened, the Calibrators are stable six months at 2-8°C.

### 6.2. Preparation of the Sample

The determination of Progesterone can be performed in human plasma as well as in serum.

Store reagent at -20°C if the determination is not performed on the same day of the sample connection. Avoid repetitive freezing and thawing of samples.

The Control is ready for use.

### 6.3. Preparation of the Wash Solution

Dilute the contents of each vial of the "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C. In concentrated wash solution it

is possible to observe the presence of crystals. In this case mix at room temperature until complete dissolution of crystals is observed. For greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer the crystals completely, then mix until the crystals are completely dissolved.

#### 6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>4</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Calibrator C <sub>0</sub> -C <sub>4</sub>	20 µL		
Sample/ Control		20 µL	
Conjugate	200 µL	200 µL	
Incubate at 37°C for 1 hour. Remove the content from each well; wash the wells 3 times with 300 µL of diluted Wash Solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature 22±28°C for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake gently the microplate. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Progesterone for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 8. RESULTS

### 8.1. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C<sub>0</sub>-C<sub>4</sub>) and of each sample.

### 8.2. Calibration curve

Plot the mean value of absorbance (Em) of the Calibrators (C<sub>0</sub>-C<sub>4</sub>) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

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Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of

expected values based on the indigenous population where the laboratory works.

Ed. 01/2015

DCM006-9

**DiaMetra S.r.l. Headquarter:** Via Calabria 15  
20090 SEGRATE (MI) Italy  
Tel. +39-02-2139184  
Fax +39-02-2133354  
**Manufactory:** Via Pozzuolo 14,  
06038 SPELLO (PG) Italy  
Tel. +39-0742-24851  
Fax +39-0742-316197  
E-mail: [info@diametra.com](mailto:info@diametra.com)

## 10. PERFORMANCE AND CHARACTERISTICS

### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (20x) of two different control sera in one assay. The within assay variability is 4%.

#### 10.1.2. Inter Assay Variation

Between run variation was determined by replicate measurements (10x) of three different control sera in different lots of kit. The between assay variability is 9.3%.

### 10.2. Accuracy

The recovery of 1.0 - 2.0 - 4.0 - 8.0 ng/mL of Progesterone added to sample gave an average value ( $\pm$ SD) of 100.88%  $\pm$  8.29% with reference to the original concentrations.

### 10.3. Sensitivity

The lowest detectable concentration of Progesterone that can be distinguished from the Calibrator 0 is 0.05 ng/mL at the 95% confidence limit.

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### 10.5. Correlation with RIA

Diametra Progesterone ELISA kit was compared to Adaltis Progesterone EIAGen kit commercially available. 31 serum samples were analysed according in both test systems.

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## BIBLIOGRAFIA

1. Kosasa T.S., "Measurement of Human Luteinizing Hormone." *Journal of Reproductive Medicine*, **26** (1981) pg. 201-6.
2. Danzer H., Braunstein G.D., et al., "Maternal Serum Human Chorionic Gonadotropic Concentrations and Fetal Sex Predictions." *Fertility and Sterility*, **34** (1980) pg. 336-40.
3. Braunstein G.D., et al., "Serum Human Luteinizing Hormone Levels through Normal Pregnancy", *American Journal of Obstetrics and Gynecology*, **126** (1976) pg. 678-81.
4. Goldstein D.P., and Kosasa T.S., "The Subunit Radioimmunoassay for LH Clinical Application." *Gynecology*, **6** (1975) pg. 145-84.
5. Batzer F., "Hormonal Evaluation of Early Pregnancy", *Fertility and Sterility*, **34** (1980) pg. 1-12.
6. Braunstein, G.D., et al., "First-Trimester Luteinizing Hormone Measurements as an Aid to the Diagnosis of Early Pregnancy Disorders", *American Journal of Obstetrics and Gynecology*, **131** (1978) pg. 25-32.
7. Lenton E., Neal L. and Sulaiman R., "Plasma Concentrations of Human Gonadotropin from the time of Implantation until the Second Week of Pregnancy", *Fertility and Sterility*, **37** (1982) pg. 773-78.

Ed. 01/2015

DCM009-12

**DiaMetra S.r.l. Headquarter:** Via Calabria 15  
20090 SEGRATE (MI) Italy  
Tel. +39-02-2139184  
Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14,  
06038 SPELLO (PG) Italy  
Tel. +39-0742-24851  
Fax +39-0742-316197  
E-mail: [info@diametra.com](mailto:info@diametra.com)

DCM009-12  
Ed. 01/2015

## LH ELISA

for routine analysis

Direct immunoenzymatic determination of the luteinizing hormone (LH) in human serum or plasma.

IVD

LOT

See external label

$\Sigma = 96$  tests

REF DKO009

### INTENDED USE

DiaMetra LH ELISA kit is a solid phase enzyme immunoassay for the quantitative determination of the luteinizing hormone (LH) concentration in human serum or plasma.

LH ELISA kit is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

Luteinizing hormone (LH) is a glycoprotein consisting of two subunits with a molecular mass of 30,000 daltons. The  $\alpha$ -subunit is similar to other pituitary hormones [follicle stimulating hormone (FSH), thyroid stimulating hormone (TSH) and chorionic gonadotropin (HCG)] while the  $\beta$ -subunit is unique. The  $\beta$ -subunit confers the biological activity to the molecule.


The  $\alpha$ -subunit consists of 89 amino acid residues while the  $\beta$ -subunit contains 129 amino acids. The carbohydrate content is between 15% and 30%. The clinical usefulness of the measurement of luteinizing hormone (LH) in ascertaining the homeostasis of fertility regulation via the hypothalamic - pituitary - gonadal axis has been well established (1,2). In addition, the advent of in vitro fertilization (IVF) technology to overcome infertility associated problems has provided the impetus for rapid improvement in LH assay methodology from the technically demanding bioassay (3) to the procedurally simple and rapid immunoenzymometric assays.

### 2. PRINCIPLE OF THE TEST

Immunoenzymometric assay.

In this method, the LH calibrators, the patient specimens and/or the controls (containing the native antigen) are first added to streptavidin coated wells. Then monoclonal biotinylated and enzyme labeled antibodies are added and the reactants mixed: these antibodies have high affinity and specificity and are directed against distinct and different epitopes of LH. A reaction between the various LH antibodies and native LH occurs in the microwells without competition or steric hindrance forming a soluble sandwich complex.

The interaction is illustrated in the following equation:

**BtAb(m)** = biotinylated monoclonal antibody (excess quantity)

**AgLH** = native antigen (variable quantity)


**EnzAb** = enzyme labeled antibody (excess quantity)

**EnzAb-AgLH-BtAb(m)** = antigen-antibodies sandwich complex

**ka** = rate constant of association

**k-a** = rate constant of dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

**Streptavidin C.W.** = streptavidin immobilized on well.

**Immobilized complex** = antibodies-antigen sandwich bound.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by a washing step. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well quantitated by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

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##### 2. Control (1 vial, 1 mL)

Control concentration is indicated on the Certificate of Analysis **REF DCE045/0903-0**

##### 3. Conjugate (1 vial, 12 mL)

Antibodies anti LH conjugated with Horseradish peroxidase (HRP) and anti LH biotinilated **REF DCE002/0902-0**

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#### 3.1. Reagents necessary not supplied

Distilled water.

#### 3.2. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

#### Note

Store all reagents at 20±8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.

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- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of LH from 5 mIU/mL to 200 mIU/mL.

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## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are ready to use, are calibrated against the international reference WHO 1st IRP 68/40 and have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
mIU/mL	0	5	25	50	100	200

For sample with concentration over 200 mIU/mL dilute the sample with the Calibrator 0.

Once opened, the Calibrators are stable 6 months at 2÷8°C.

### 6.2. Preparation of the Wash solution

Dilute the contents of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2÷8°C.

### 6.3. Preparation of the sample

The specimens shall be serum or plasma; the usual precautions in the collection of venipuncture samples should be observed.

For accurate comparison to established normal values, a fasting morning serum or plasma sample should be obtained.

To obtain serum samples, the blood should be collected in a venipuncture tube without additives or anti-coagulants; allow the blood to clot; centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2÷8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, they may be stored at -20°C for up to 30 days. Avoid repetitive freezing and thawing.

When assayed in duplicate, 0,040 mL of the specimen is required.

### 6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Sample/Control		20 µL	
Calibrator C <sub>0</sub> -C <sub>5</sub>	20 µL		
Conjugate	100 µL	100 µL	
Incubate at room temperature (22÷28°C) for 1 hour. Remove the contents from each well and wash the wells 3 times with 300 µL of diluted Wash Solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22÷28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at levels in the low, medium and high ranges of the dose response curve for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

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### 5.2. Calibration Curve

Plot the mean value of absorbance (E<sub>m</sub>) of the calibrators (C<sub>0</sub>-C<sub>4</sub>) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

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## 9. EXPECTED VALUES

Reference ranges are hereby reported:

	LH (mIU/mL)
<b>MALE:</b>	0.7 – 7.4
<b>FEMALE:</b>	
<b>Follicular phase</b>	0.5 – 10.5
<b>Ovulation phase</b>	18.4 – 61.2
<b>Lutheal phase</b>	0.5 – 10.5
<b>Menopausal</b>	8.2 – 40.8

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## 10. LIMITATIONS OF THE PROCEDURE

LH is suppressed by estrogen but in woman taking oral contraceptives the level may be low or normal. Excessive dieting and weight loss may lead to low gonadotropin concentrations. Luteinizing hormone is dependent upon diverse factors other than pituitary homeostasis. Thus, the determination alone is not sufficient to assess clinical status.

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Precision

#### 11.1.1. Intra Assay Variation

Within run variation was determined by replicate (20x) the measurement of three different control sera in one assay. The within assay variability is 9.21%.

#### 11.1.2. Inter Assay Variation

Between run variation was determined by replicate (10x) the measurements of three different control sera in different lots. The between assay variability is 7.91%.

### 11.2. Correlation

Diametra LH ELISA kit was compared to a commercially available LH kit. 36 serum samples were tested.

The regression curve is:

$$\text{Diametra} = 0.91 * (\text{commercial kit}) + 0.05$$

$$R \text{ squared} = 0.971$$

The new Diametra LH kit was compared to the old Diametra LH kit. 36 serum samples were tested.

The regression curve is:

$$(\text{Diametra new}) = 1.08 * (\text{Diametra old}) - 1.22$$

$$R \text{ squared} = 0.981$$

### 11.3. Sensitivity

The lowest detectable concentration of LH that can be distinguished from the Calibrator 0 is 0.22 mIU/mL.

## 11.4. Specificity

The cross-reactivity of the Diametra LH kit to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of Luteinizing Hormone needed to produce the same absorbance.

Substance tested	Cross Reactivity
<b>LH</b>	100 %
<b>HCG</b>	0.007 %
<b>Prolactin</b>	N.D.
<b>FSH</b>	N.D.
<b>TSH</b>	N.D.

## 11.5. Accuracy

The recovery test performed on three different samples, enriched with 5.63 - 11.25 - 22.5 - 45 - 90 mIU/mL of LH, gave a average value ( $\pm$ SD) of 97.17%  $\pm$  4.00%.

In the dilution test three different samples were diluted 2, 4, 8 and 16 times with Calibrator 0; the average value ( $\pm$ SD) obtained is 99.13%  $\pm$  7.37%.

## 11.6. Hook Effect

Diametra LH assay shows no Hook Effect up to 400 mIU/mL.

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## BIBLIOGRAPHY

- 1 Kosasa T.S., "Measurement of Human Luteinizing Hormone." *Journal of Reproductive Medicine*, **26** (1981) pg. 201-6.
- 2 Danzer H., Braunstein G.D., et al., "Maternal Serum Human Chorionic Gonadotropic Concentrations and Fetal Sex Predictions." *Fertility and Sterility*, **34** (1980) pg. 336-40.
- 3 Braunstein G.D., et al., "Serum Human Luteinizing Hormone Levels through Normal Pregnancy", *American Journal of Obstetrics and Gynecology*, **126** (1976) pg. 678-81.
- 4 Goldstein D.P., and Kosasa T.S., "The Subunit Radioimmunoassay for LH Clinical Application." *Gynecology*, **6** (1975) pg. 145-84.
- 5 Batzer F., "Hormonal Evaluation of Early Pregnancy", *Fertility and Sterility*, **34** (1980) pg. 1-12.
- 6 Braunstein, G.D., et al., "First-Trimester Luteinizing Hormone Measurements as an Aid to the Diagnosis of Early Pregnancy Disorders", *American Journal of Obstetrics and Gynecology*, **131** (1978) pg. 25-32.
- 7 Lenton E., Neal L. and Sulaiman R., "Plasma Concentrations of Human Gonadotropin from the time of Implantation until the Second Week of Pregnancy", *Fertility and Sterility*, **37** (1982) pg. 773-78.

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1. Kosasa T.S., "Measurement of Human Luteinizing Hormone." *Journal of Reproductive Medicine*, **26** (1981) pg. 201-6.
2. Danzer H., Braunstein G.D., et al., "Maternal Serum Human Chorionic Gonadotropic Concentrations and Fetal Sex Predictions." *Fertility and Sterility*, **34** (1980) pg. 336-40.
3. Braunstein G.D., et al., "Serum Human Luteinizing Hormone Levels through Normal Pregnancy", *American Journal of Obstetrics and Gynecology*, **126** (1976) pg. 678-81.
4. Goldstein D.P., and Kosasa T.S., "The Subunit Radioimmunoassay for LH Clinical Application." *Gynecology*, **6** (1975) pg. 145-84.
5. Batzer F., "Hormonal Evaluation of Early Pregnancy", *Fertility and Sterility*, **34** (1980) pg. 1-12.

6. Braunstein, G.D., et al., "First-Trimester Luteinizing Hormone Measurements as an Aid to the Diagnosis of Early Pregnancy Disorders", *American Journal of Obstetrics and Gynecology*, **131** (1978) pg. 25-32.
7. Lenton E., Neal L. and Sulaiman R., "Plasma Concentrations of Human Gonadotropin from the time of Implantation until the Second Week of Pregnancy ", *Fertility and Sterility*, **37** (1982) pg. 773-78.

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## BIBLIOGRAFIA

- Odell, W.D., Parlow, A.F., et al, *J Clin Invest*, 47, 2551 (1981).
- Saxema, B.B., Demura, H.M., et al, *J Clin Endocrinol Metab.*, 28, 59 (1968).
- Wennink JM, Delemarre-van de Waal HA, Schoemaker R, Schoemaker H, Schoemaker J. 1990 Luteinizing hormone and follicle stimulating hormone secretion patterns in girls throughout puberty measured using highly sensitive immunoradiometric assays. *Clin Endocrinol (Oxf)*. 33:333–344.
- Winter JS, Faiman C. 1973 The development of cyclic pituitarygonadal function in adolescent females. *J Clin Endocrinol Metab.* 37:714–718.
- Simoni M, Gromoll J, Nieschlag E 1997 The follicle stimulating hormone receptor: biochemistry, molecular biology, physiology and pathophysiology. *Endocr Rev* 18:739–773.
- Vitt UA, Kloosterboer HJ, Rose UM, Mulders JW, Kiesel PS, Bete S, Nayudu PL 1998 Isoforms of human recombinant folliclestimulating hormone: comparison of effects on murine follicle development *in vitro*. *Biol Reprod* 59:854–861.
- Layman LC, Lee EJ, Peak DB, Namnoum AB, Vu KV, van Lingen BL, Gray MR, McDonough PG, Reindollar RH, Jameson JL 1997 Delayed puberty and hypogonadism caused by mutations in the folliclestimulating hormone  $\beta$  subunit gene. *N Engl J Med* 337:607–611.
- Robertson DR 1991 Circulating half-lives of follicle stimulating hormones and luteinizing hormone in pituitary extracts and isoform fractions of ovariectomized and intact ewes. *Endocrinology* 129:1805– 1813.
- Wide L 1981 Electrophoretic and gel chromatographic analyses of follicle stimulating hormone in human serum. *Ups J Med Sci* 86:249–258.
- Berger P, Bidart JM, Delves PS, Dirrhofer S, Hoermann R, Isaacs N, Jackson A, Klonisch T, Laphorn A, Lund T, Mann K, Roitt I, SchwarzS, Wick G 1996 Immunochemical mapping of gonadotropins. *Mol Cell Endocrinol* 125:33–43.

Ed. 01/2015

DCM010-13

**DiaMetra S.r.l. Headquater:** Via Calabria 15

20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14,

06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: [info@diametra.com](mailto:info@diametra.com)

DCM010-13  
Ed. 01/2015

## FSH ELISA

for routine analysis

Direct immunoenzymatic determination of the Follicle-Stimulating Hormone (FSH) in human serum or plasma.

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO010

### INTENDED USE

Diametra FSH ELISA kit is a direct solid phase immunoassay for the quantitative determination of Follicle-Stimulating Hormone (FSH) in human serum or plasma.

FSH ELISA kit is intended for laboratory use only.

### 1. SUMMARY AND EXPLANATION OF THE TEST

Follicle Stimulating hormone (FSH) is a glycoprotein consisting of two subunits with an approximate molecular mass of 35,500 daltons. The  $\alpha$ - subunit is similar to other pituitary hormones [luteinizing stimulating hormone (LH), thyroid stimulating hormone (TSH) and chorionic gonadotropin (hCG)] while the  $\beta$ -subunit is unique. The  $\beta$ - subunit confers the biological activity to the molecule. Stimulation by gonadotropin-releasing hormone (GnRH) causes release of FSH, as well as LH, from the pituitary and is transported by the blood to their sites of action, the testes or ovary.

In men, FSH acts on the Sertoli cells of the testis, stimulating the synthesis of inhibin, which appears to specifically inhibit further FSH secretion, and androgen-binding protein. Thus, it indirectly supports spermatogenesis.

In women, FSH acts on the granulosa cells of the ovary, stimulating steroidogenesis. All ovulatory menstrual cycles have a characteristic pattern of FSH, as well as LH, secretion. The menstrual cycle is divided into a follicular phase and a luteal phase by the midcycle surge of the gonadotropins (LH and FSH). As the follicular phase progresses, FSH concentration decreases. Near the time ovulation occur, about midcycle, FSH peaks (lesser in magnitude than LH) to its highest level.


The clinical usefulness of the measurement of Follicle Stimulating hormone (FSH) in ascertaining the homeostasis of fertility regulation via the hypothalamic - pituitary - gonadal axis has been well established (1,2).

### 2. PRINCIPLE OF THE TEST

The essential reagents required for an immunoenzymatic assay include high affinity and specificity antibodies (enzyme-linked and immobilised) with different and distinct epitope recognition, in excess, and native antigen.

In this procedure the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on the well and exogenously added biotinylated monoclonal anti FSH antibody.

Upon mixing monoclonal biotinylated antibody, the enzyme labeled antibody and a serum containing the native antigen, reaction results between the native antigen and the antibodies without competition or steric hindrance to form a soluble sandwich complex. The interaction is illustrated by the following equation:

$\text{B}^{\text{tn}}\text{Ab}_{(m)}$  = biotinylated monoclonal antibody (Excess quantity)

$\text{Ag}_{\text{FSH}}$  = native FSH antigen (variable quantity)

$\text{EnzAb}_{(p)}$  = enzyme labeled polyclonal antibody (Excess quantity)


$\text{EnzAb}_{(p)} - \text{Ag}_{\text{FSH}} - \text{B}^{\text{tn}}\text{Ab}_{(m)}$  = antigen-antibodies sandwich complex

$\text{K}_a$  = rate constant of association

$\text{K}_{-a}$  = rate constant of disassociation

Simultaneously the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody.

This interaction is illustrated below:

$\text{Streptavidin}_{\text{cw}}$  = streptavidin immobilized on well

Immobilized complex = antibodies-antigen sandwich bound.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by a washing step.

The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration.

By using several different serum references of known antigen values, a dose response curve can be

generated from which the antigen concentration of an unknown can be ascertained.

### 3. REAGENTS, MATERIALS AND INSTRUMENTATION

#### 3.1 Reagents and Materials Supplied in the Kit

##### 1. Calibrators (6 vials, 1 mL each)

CAL0	<b>REF</b> DCE002/1006-0
CAL1	<b>REF</b> DCE002/1007-0
CAL2	<b>REF</b> DCE002/1008-0
CAL3	<b>REF</b> DCE002/1009-0
CAL4	<b>REF</b> DCE002/1010-0
CAL5	<b>REF</b> DCE002/1011-0

##### 2. Control (1 vial, 1 mL)

Control concentration is indicated on the Certificate of Analysis **REF** DCE045/1003-0

##### 3. Conjugate (1 vial, 12 mL)

Antibody anti FSH conjugated with Horseradish peroxidase (HRP)

Antibody anti FSH conjugated with Biotine

**REF** DCE002/1002-0

##### 4. Coated Microplate (1 breakable microplate)

Streptavidin adsorbed on microplate

**REF** DCE002/1003-0

##### 5. TMB Substrate (1 vial, 15 mL)

H<sub>2</sub>O<sub>2</sub>-TMB 0.26 g /L (avoid any skin contact)

**REF** DCE004-0

##### 6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)

**REF** DCE005-0

##### 7. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M pH 7.4 **REF** DCE054-0

#### 3.2 Reagents necessary not supplied

Distilled water.

#### 3.3 Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

#### Note

Store all reagents at 2-8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable up to expiry date of the kit.

#### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.

- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of FSH from 5 to 100 mIU/mL.

#### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are ready to use, are calibrated against the International Standard WHO 2<sup>nd</sup> IRP 78/549 and have the following concentration:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
mIU/mL	0	5	10	25	50	100

For sample with concentration over 100 mIU/mL dilute the sample 1:2 with the Calibrator 0.

Once opened, the Calibrators are stable 6 months at 2±8°C.

### 6.2. Preparation of the Wash Solution

Dilute the content of each vial of the "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2±8°C.

In concentrated wash solution is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of crystals; for greater accuracy dilute the whole content of the bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely, then mix until crystals are completely dissolved.

### 6.3. Preparation of the sample

The specimens should be human serum or plasma.

Note: it is recommended to avoid the use of EDTA-plasma, because it could interfere with the sample assay.

The usual precautions in the collection of venipuncture samples should be observed.

For accurate comparison to established normal values, a fasting morning serum sample should be obtained.

To obtain the serum, the blood should be collected in a venipuncture tube without additives or anti-coagulants; allow the blood to clot; centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2±8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, the samples may be stored at -20°C for up to 30 days. Avoid repetitive freezing and thawing.

When assayed in duplicate, 0.100 mL of the specimen is required.

### 6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the

calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Sample /Control		50 µL	
Calibrator C <sub>0</sub> -C <sub>5</sub>	50 µL		
Conjugate	100 µL	100 µL	
Incubate at room temperature (22-28°C) for 1 hours. Remove the contents from each well. Wash the wells 3 times with 300 µL of diluted Wash Solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22-28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at levels of a low, normal, and high range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

If computer controlled data reduction is used to calculate the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

## 8. RESULTS

### 8.1. Mean Absorbance

Calculate the mean of the absorbance (E<sub>m</sub>) for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>) and of each sample.

### 8.2. Calibration Curve

Plot the mean value of absorbance (E<sub>m</sub>) of the calibrators (C<sub>0</sub>-C<sub>5</sub>) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

### 8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in mIU/mL.

## 9. EXPECTED VALUES

Reference ranges are reported below:

	FSH (mIU/mL)
<b>MALE</b>	1 – 14
<b>FEMALE:</b>	
Follicular phase	3 – 12
Mid-cycle	8 – 22
Luteal phase	2 – 12
Menopausal	35 – 151

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

## 10. LIMITATIONS OF THE PROCEDURE

- **For diagnostic purposes FSH values should be used as an adjunct to other data available to the physician.** Procedural directions must be followed exactly and careful technique must be used to obtain valid results. Any modification of the procedure is likely to alter the results. FSH is dependent upon diverse factors other than pituitary homeostasis. Thus, the determination of FSH alone is not sufficient to assess clinical status.
- FSH is suppressed by estrogen but in woman taking contraceptives the level may be low or normal. Excessive dieting and weight loss may lead to low gonadotropin concentration.
- Patients specimens with abnormally high FSH levels can cause a hook effect, that is, paradoxical low absorbance results. If this is suspected, dilute the specimen 1:100 with the Calibrator 0, re-assay (multiply the result by 100). However, values as high as 2000 mIU/mL have been found to absorb greater than the absorbance of the highest calibrator.
- Patients receiving preparations of mouse monoclonal antibodies for diagnosis of therapy may contain human anti-mouse antibodies (HAMA) and may show either falsely elevated or depressed values when assayed.

## 11. PERFORMANCE AND CHARACTERISTICS

### 11.1. Precision

#### 11.1.1. Intra-Assay

Within-run precision was determined by replicate (20x) the measurement of three different control sera in one assay. The within-assay variability is  $\leq 9.7\%$ .

#### 11.1.2. Inter-Assay

Between-run precision was determined by replicate (10x) the measurement of three different control sera in different lots of kits. The between-assay variability is  $\leq 8.0\%$ .

### 11.2. Correlation

Diametra FSH kit (Y) was compared to a commercially available FSH kit (X). 35 serum samples have been tested. The linear regression curve was:

$$Y = 0.91 \cdot X + 3.41$$

$$r^2 = 0.969$$

The new Diametra FSH kit (Y) was compared to the previous Diametra FSH kit (X): 60 serum samples have been tested. The linear regression curve was:

$$Y = 1.01 \cdot X - 0.44$$

$$r^2 = 0.953$$

### 11.3. Sensitivity

The minimal detectable concentration of FSH that can be distinguished from the Calibrator 0 is 0.17 mIU/mL.

### 11.4. Specificity

The cross-reactivity of the Diametra FSH ELISA method to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of Follicle Stimulating Hormone needed to produce the same OD.

Substance	Concentration	Cross Reactivity
<b>FSH</b>	---	100,0%
<b>HCG</b>	5000 mIU/mL	N.D.
<b>LH</b>	400 mIU/mL	N.D.
<b>Prolactin</b>	2000 ng/mL	N.D.
<b>TSH</b>	1000 mIU/L	N.D.

### 11.5. Accuracy

Three different samples were spiked with 3, 6, 12 and 24 mIU/mL of FSH; the recovery test showed an average value ( $\pm$  SD) of  $97.83\% \pm 8.52\%$ .

Three different samples were diluted 2, 4 and 8 times with the Calibrator 0; the dilution test showed an average value ( $\pm$  SD) of  $105.54\% \pm 5.36\%$ .

## 12. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

## BIBLIOGRAPHY

- Odell, W.D., Parlow, A.F., et al, *J Clin Invest*, 47, 2551(1981).
- Saxema, B.B., Demura, H.M., et al, *J Clin Endocrinol Metab.*, 28, 591 (1968).
- Wennink JM, Delemarre-van de Waal HA, Schoemaker R, Schoemaker H, Schoemaker J. 1990 Luteinizing hormone and follicle stimulating hormone secretion patterns in girls throughout puberty measured using highly sensitive immunoradiometric assays. *Clin Endocrinol (Oxf)* 33:333–344.
- Winter JS, Faiman C. 1973 The development of cyclic pituitarygonadal function in adolescent females. *J Clin Endocrinol Metab.* 37:714–718.
- Simoni M, Gromoll J, Nieschlag E 1997 The follicle stimulating hormone receptor: biochemistry, molecular biology, physiology and pathophysiology. *Endocr Rev* 18:739–773.
- Vitt UA, Kloosterboer HJ, Rose UM, Mulders JW, Kiesel PS, Bete S, Nayudu PL 1998 Isoforms of human recombinant folliclestimulating hormone: comparison of effects on murine follicle development *in vitro*. *Biol Reprod* 59:854–861.
- Layman LC, Lee EJ, Peak DB, Namnoum AB, Vu KV, van Lingen BL, Gray MR, McDonough PG, Reindollar RH, Jameson JL 1997 Delayed puberty and hypogonadism caused by mutations in the follicle-stimulating hormone  $\beta$  subunit gene. *N Engl J Med* 337:607–611.
- Robertson DR 1991 Circulating half-lives of follicle stimulating hormones and luteinizing hormone in pituitary extracts and isoform fractions of ovariectomized and intact ewes. *Endocrinology* 129:1805–1813.
- Wide L 1981 Electrophoretic and gel chromatographic analyses of follicle stimulating hormone in human serum. *Ups J Med Sci* 86:249–258.
- Berger P, Bidart JM, Delves PS, Dirnhofer S, Hoermann R, Isaacs N Jackson A, Klonisch T, Laphorn A, Lund T, Mann K, Roitt I, Schwarz S, Wick G 1996 Immunochemical mapping of gonadotropins. *Mol Cell Endocrinol* 125:33–43.

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- Odell, W.D., Parlow, A.F., et al, *J Clin Invest*, 47, 2551 (1981).
- Saxema, B.B., Demura, H.M., et al, *J Clin Endocrinol Metab.*, 28, 59 (1968).
- Wennink JM, Delemarre-van de Waal HA, Schoemaker R, Schoemaker H, Schoemaker J. 1990 Luteinizing hormone and follicle stimulating hormone secretion patterns in girls throughout puberty measured using highly sensitive immunoradiometric assays. *Clin Endocrinol (Oxf)*. 33:333–344.
- Winter JS, Faiman C. 1973 The development of cyclic pituitarygonadal function in adolescent females. *J Clin Endocrinol Metab*. 37:714–718.
- Simoni M, Gromoll J, Nieschlag E 1997 The follicle stimulating hormone receptor: biochemistry, molecular biology, physiology and pathophysiology. *Endocr Rev* 18:739–773.
- Vitt UA, Kloosterboer HJ, Rose UM, Mulders JW, Kiesel PS, Bete S, Nayudu PL 1998 Isoforms of human recombinant folliclestimulating hormone: comparison of effects on murine follicle development *in vitro*. *Biol Reprod* 59:854–861.
- Layman LC, Lee EJ, Peak DB, Namnoum AB, Vu KV, van Lingen BL, Gray MR, McDonough PG, Reindollar RH, Jameson JL 1997 Delayed puberty and hypogonadism caused by mutations in the folliclestimulating hormone  $\beta$  subunit gene. *N Engl J Med* 337:607–611.
- Robertson DR 1991 Circulating half-lives of follicle stimulating hormones and luteinizing hormone in pituitary extracts and isoform fractions of

ovariectomized and intact ewes. *Endocrinology* 129:1805– 1813.

- Wide L 1981 Electrophoretic and gel chromatographic analyses of follicle stimulating hormone in human serum. *Ups J Med Sci* 86:249–258.
- Berger P, Bidart JM, Delves PS, Dirnhofer S, Hoermann R, Isaacs N, Jackson A, Klonisch T, Laphorn A, Lund T, Mann K, Roitt I, Schwarz S, Wick G 1996 Immunochemical mapping of gonadotropins. *Mol Cell Endocrinol* 125:33–43.

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# PROLACTIN ELISA

for routine analysis

Direct immunoenzymatic determination of Prolactin in human serum.

IVD

LOT

See external label

Σ = 96 tests

REF DKO011

## INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of Prolactin concentration in human serum.

Prolactin ELISA kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

Prolactin is a polypeptide hormone synthesized and secreted by the Adenohypophysis (anterior Pituitary gland) and the placenta. It is also produced in other tissues including the breast and the decidua. Pituitary prolactin secretion is regulated by neuroendocrine neurons in the hypothalamus, most importantly by neurosecretory dopamine neurons of the arcuate nucleus, which inhibit prolactin secretion.

Prolactin is present in several body fluids, including blood plasma, amniotic fluid, milk, mucosal secretions and cerebrospinal fluid.

Prolactin has many effects, the most important of which is to stimulate the mammary glands to produce milk (lactation). Other possible functions of prolactin include the surfactant synthesis of the fetal lungs at the end of the pregnancy and immune tolerance of the foetus by the maternal organism during pregnancy.

Prolactin may also have inhibitory effects on gonadal function when present in high concentrations.

There is a diurnal cycle in prolactin secretion.

During pregnancy, high circulating concentrations of estrogen promote prolactin production. The resulting high levels of prolactin secretion cause maturation of the mammary glands, preparing them for lactation. After childbirth, prolactin levels fall as the internal stimulus for them is removed.

High prolactin levels also tend to suppress the ovulatory cycle by inhibiting the secretion of both FSH and GnRH.

Prolactin levels may be checked as part of a sex hormone workup, as elevated prolactin secretion can suppress the secretion of FSH and GnRH, leading to hypogonadism, and sometimes causing erectile dysfunction in men.

Elevations in plasma prolactin concentrations occur during ovulation, pregnancy, nursing and stress. Abnormal elevations in plasma prolactin levels (hyperprolactinemia) can occur as a result of pituitary adenomas, other anatomic and traumatic

abnormalities, in response to certain pharmacologic agents and in hypothyroidism.

Hypoprolactinemia (low prolactin levels) are observed in cases of hypopituitarism.


## 2. PRINCIPLE

The essential reagents required for an immunoenzymatic assay include high affinity and specificity antibodies (enzyme and immobilised) with different and distinct epitope recognition, in excess, and native antigen.

In this procedure the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on the well and exogenously added biotinylated monoclonal anti-PRL antibody.

Upon mixing monoclonal biotinylated antibody, the enzyme labeled antibody and a serum containing the native antigen, a reaction results between the native antigen and the antibodies without competition or steric hindrance to form a soluble sandwich complex.

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$\text{B}^{\text{tn}}\text{Ab}_{(m)}$  = Biotinylated Monoclonal Antibody (Excess quantity)

$\text{Ag}_{\text{PRL}}$  = Native PRL antigen (variable quantity)

$\text{EnzAb}_{(p)}$  = Enzyme labeled polyclonal antibody (Excess quantity)


$\text{EnzAb}_{(p)} - \text{Ag}_{\text{PRL}} - \text{B}^{\text{tn}}\text{Ab}_{(m)}$  = Antigen-Antibodies Sandwich complex

$\text{K}_a$  = Rate constant of association

$\text{K}_a$  = Rate constant of disassociation

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Streptavidin<sub>cw</sub> = Streptavidin immobilized on well  
Immobilized complex = Antibodies-Antigen sandwich bound.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. By using several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

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##### 2. Control (1 vial, 1 mL)

Control Concentration is indicated on the Certificate of Analysis **REF** **DCE045/1103-0**

##### 3. Conjugate (1 vial, 12 mL)

Antibodies Anti Prolactin conjugated with Horseradish peroxidase (HRP) and Anti Prolactin biotinylated **REF** **DCE002/1102-0**

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Streptavidin adsorbed on microplate **REF** **DCE002/1103-0**

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### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of Prolactin from 5.0 to 100.0 ng/mL.

### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.

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## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are ready to use, are calibrated against WHO 3<sup>rd</sup> IS 84/500 and have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
ng/mL	0	5	10	25	50	100

The Calibrators are stable until the expiry date printed on the label. Once opened, the Calibrators are stable 6 months at 2-8°C.

### 6.2. Preparation of Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (10X) with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C. In concentrated wash solution is possible to observe the presence of crystals; in this case mix at room temperature until the complete dissolution of crystals; for greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care to transfer completely the crystals, then mix until crystals are completely dissolved.

### 6.3. Preparation of the Sample

Prolactin determination can be done on human serum. **Do not use plasma for this assay** (plasma samples can lead to false results).

Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.

To obtain the serum, collect blood by venipuncture, allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred.

Specimen can be stored at 2÷8°C for at short time (max two days). For longer storage the specimen should be frozen at -20°C. Avoid repeated freezing and thawing. Thawed samples should be inverted several times prior to testing.

For sample with concentration over 100 ng/mL dilute the sample 1:2 with Calibrator 0.

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- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Sample/Control		50 µL	
Calibrator C <sub>0</sub> -C <sub>5</sub>	50 µL		
Conjugate	100 µL	100 µL	
Incubate at room temperature (22÷28°C) for 1 hours. Remove the contents from each well; wash the wells 3 times with 300 µL of diluted Wash Solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22÷28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Prolactin for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

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### 8.2. Calibration curve

Plot the values of absorbance (Em) of the calibrators (C<sub>0</sub>-C<sub>5</sub>) against concentration. Draw the best-fit curve through the plotted points (Es: Four Parameter Logistic or Sigmoide).

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## 9. REFERENCE VALUES

Each laboratory must establish its own normal ranges based on patient population.

The serum Prolactin values are comprised in the following intervals:

Sample		Range ng/mL
Male		1.8 - 17.0
Female:	menstrual cycle	1.2 - 19.5
	menopause	1.5 - 18.5

Some of the female population tested in this group were probably using oral contraceptives, which may affect results.

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### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate (20x) the measurement of three different control sera in one assay.

Sample	N	X	σ	C.V.
Level 1	20	5.33	0.15	2.78%
Level 2	20	18.21	0.73	4.03%
Level 3	20	37.20	1.38	3.71%

#### 10.1.2. Inter Assay Variation

Between run variation was determined by replicate (10x) the measurement of three different control sera with kits of different lots.

Sample	N	X	σ	C.V.
Level 1	10	5.46	0.30	5.49%
Level 2	10	17.72	0.91	5.16%
Level 3	10	36.29	1.67	4.60%

### 10.2. Sensitivity

The lowest detectable concentration of Prolactin that can be distinguished from the Calibrator 0 is 0.12 ng/mL at the 95 % confidence limit.

### 10.3. Accuracy

The recovery of 3,13 - 6,25 - 12,50 - 25,00 - 50,00 ng/mL of Prolactin added to sample gave an average value (±SD) of 102.52% ± 9.75% with reference to the original concentrations.

The dilution test performed on three sera diluted 2 - 4 - 8 - 16 times gave an average value (±SD) of 102.19% ± 9.80%.

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### 10.5. Correlation

Diametra Prolactin ELISA kit (y) was compared to another commercially available Prolactin assay (x). Serum samples of 37 subjects were analysed.

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$$y = 1.01 x + 1.94$$

$$r^2 = 0.957$$

The new Diametra Prolactin ELISA kit (y) was compared to the old Diametra Prolactin ELISA kit (x). Serum samples of 37 subjects were analysed.

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$$r^2 = 0.937$$

### 10.6. Hook Effect

Diametra Prolactin assay shows no Hook effect up to 200 ng/mL.

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## AFP ELISA

for routine analysis

Direct immunoenzymatic determination of AFP in human serum or plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO012

### INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of AFP concentration in human serum or plasma.

AFP ELISA kit is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

Alpha Fetoprotein (AFP) is a 68 kDa glycoprotein, which is normally only produced in the fetus during its development. It is a normally produced by the liver and yolk sac of the fetus.

AFP levels decrease soon after birth and probably has no function in normal adults. It binds the hormone estradiol to keep it from affecting the fetal brain. Its measurement during pregnancy has been useful to detect certain abnormalities - specifically, if high levels of AFP are found in amniotic fluid, it can indicate a developmental defect in the baby.

In some patients who are not pregnant a tumor can produce AFP, thus it can be used as a tumour marker. AFP is the main tumour marker (along with HCG) to diagnose testicular cancer and its values over time can have significant effect on the treatment plan.

Like all tumour markers, the detection of AFP by itself is not diagnostic of anything, although if it is detected it is certainly advisable to rule out the diseases could cause levels to rise. The primary reason tumor markers are used are to measure the success of a treatment (e.g. chemotherapy), if levels of AFP are going down, it is an indication that a disease is improving. New research exhibits that an isoform of AFP which binds *Lens culinaris* agglutinin (AFP-L3) can be particularly useful in early identification of aggressive tumors associated with hepatocellular carcinoma (HCC).

### 2. PRINCIPLE

AFP ELISA test is based on simultaneous binding of human AFP to two monoclonal antibodies, one immobilized on microwell plates, the other conjugated with horseradish peroxidase (HRP).

After incubation, the separation bound-free is obtained with a simple solid-phase washing.

Then the enzyme HRP in the bound-fraction reacts with the Substrate (H<sub>2</sub>O<sub>2</sub>) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H<sub>2</sub>SO<sub>4</sub>) is added.

The colour intensity is proportional to the AFP concentration in the sample.

The AFP concentration in the sample is calculated through a calibration curve.

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CAL4 REF DCE002/1210-0

2. Control (1 vial, 1 mL)

Concentration of Control is indicated on the Certificate of Analysis REF DCE045/1203-0

3. Incubation Buffer (1 vial, 50 mL)

Phosphate buffer 50 mM pH 7.4; BSA 1 g/L

REF DCE001/1401-0

4. Conjugate (1 vial, 1 mL)

Monoclonal anti AFP antibody conjugated with Horseradish peroxidase (HRP) REF DCE002/1202-0

5. Coated Microplate (1 breakable microplate)

Monoclonal anti AFP antibody adsorbed on microplate REF DCE002/1203-0

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## Note

Store all reagents at 20±8°C in the dark.

Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until expiry date of the kit.

## 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagents should be handled in the same manner as potentially infectious material.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of AFP from 5 ng/mL to 200 ng/mL.

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## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>4</sub>)

The Calibrators are ready to use, are calibrated against the WHO 1<sup>st</sup> IS 72/225 and have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>
ng/mL	0	5	20	80	200

The Calibrators are stable until the expiry date printed on the label. Once opened, the calibrators are stable six months at 2-8°C

### 6.2. Preparation of Diluted Conjugate

Prepare immediately before use.

Add 10 µL of Conjugate (reagent 4) to 1 mL of Incubation Buffer (reagent 3). The quantity of diluted conjugate is proportional at the number of the tests. Mix gently for 5 minutes, with a rotating mixer. Stable for 3 hours at room temperature (22±28°C).

### 6.3. Preparation of the Sample

AFP determination should be done in human serum or plasma.

Specimen can be stored at 2±8°C for at short time (max two days). For longer storage the specimen should be frozen at -20°C. Avoid repeated freezing and thawing.

For sample with concentration over 200 ng/mL dilute the sample with Incubation buffer.

The Control is ready to use.

#### 6.4. Preparation of Wash Solution

Dilute the content of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

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Reagent	Calibrator	Sample/ Control	Blank
Calibrator C <sub>0</sub> -C <sub>4</sub>	25 µL		
Sample/ Control		25 µL	
Diluted Conjugate	50 µL	50 µL	
Incubate at room temperature (22-28°C) for 1 hour. Remove the contents from each well; wash the wells 3 times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22-28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

#### 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of AFP for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be

employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

#### 8. RESULTS

##### 8.1. Mean Absorbance

Calculate the mean of the absorbance (E<sub>m</sub>) for each point of the calibration curve (C<sub>0</sub>-C<sub>4</sub>) and of each sample.

##### 8.2. Calibration curve

Plot the values of absorbance (E<sub>m</sub>) of the calibrators (C<sub>0</sub>-C<sub>4</sub>) against concentration. Draw the best-fit curve through the plotted points. (Es: Four Parameter Logistic).

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#### 9. REFERENCE VALUES

In a study conducted with apparently normal healthy adults, using Diametra AFP, the following results were observed:

Population	0-10 ng/mL	20 ng/mL	30 ng/mL
<b>Males</b>	82	2	1
<b>Females</b>	55	1	1

In a study conducted with nonseminomatous testicular cancer patient using AFP, the following values were observed:

Population	0-10 ng/mL	10-100 ng/mL	>100 ng/mL
<b>Males</b>	4	5	3

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### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate (16x) the measurement of three different control sera in one assay. The within assay variability is 6.4%.

#### 10.1.2. Inter Assay Variation

Between run variation was determined by replicate (16x) the measurements of three different control sera in different lots. The between assay variability is 6.5%.

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### 10.3. Accuracy

The recovery of 12.5 – 25 – 50 ng/mL of AFP added to sample gave an average value ( $\pm$ SD) of 98.33%  $\pm$  4.43% with reference to the original concentrations.

### 10.4. Sensitivity

The lowest detectable concentration of AFP that can be distinguished from the Calibrator 0 is 0.35 ng/mL at the 95% confidence limit.

### 10.5. Correlation with RIA

Diametra AFP ELISA was compared to another commercially available AFP assay. Serum samples of 32 females and 4 males were analysed according in both test systems.

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$$r = 0.99 (r^2 = 0.98)$$

### 10.6. Hook Effect

The AFP ELISA, a competitive enzyme immunoassay, shows no Hook Effect up to 4000 ng/mL.

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Tel. +39-02-2139184

Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14  
06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: [info@diametra.com](mailto:info@diametra.com)

DCM013-11  
Ed. 01/2015

## TSH ELISA

for routine analysis


Direct immunoenzymatic determination of TSH in human serum or plasma

IVD

LOT

See external label

Σ = 96 tests

REF DKO013

### INTENDED USE

Diametra TSH ELISA is an immunoenzymatic colorimetric method for quantitative determination of thyroid-stimulating hormone (TSH, thyrotropin) concentration in human serum or plasma.

TSH ELISA kit is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

Thyroid-stimulating hormone (TSH or thyrotropin) is a glycopolypeptide hormone synthesized and secreted by the anterior pituitary gland which regulates the endocrine function of the thyroid gland.

TSH consists of two subunits, the *alpha* and the *beta* subunit. The  $\alpha$  subunit is identical to that of human chorionic gonadotropin (HCG), luteinising hormone (LH), follicle-stimulating hormone (FSH). The (*beta*) subunit is unique to TSH, and therefore determines its function.

TSH production is controlled by a Thyrotropin Releasing Hormone, (TRH), which is manufactured in the hypothalamus and transported to the pituitary gland, where it increases TSH production and release. Somatostatin is also produced by the hypothalamus, and has an opposite effect on the pituitary production of TSH, decreasing or inhibiting its release.

TSH stimulates the thyroid gland to secrete the hormones thyroxine (T4) and triiodothyronine (T3).

Release of TSH is regulated by the circulating free fraction of thyroid hormones in the blood. TSH levels are depressed when peripheral concentrations of the free fraction of thyroid hormones are high. Conversely, TSH levels are high when peripheral concentrations of thyroid hormones are low. This effect creates a regulatory negative feedback loop.

TSH levels are tested in the blood of patients suspected of suffering from excess (hyperthyroidism), or deficiency (hypothyroidism) of thyroid hormone. Generally, a normal range for TSH is between 0.3 and 3.0 mIU/mL, but the interpretation depends also on what the blood levels of thyroid hormones (T3 and T4) are.

Higher than normal levels of TSH combined with high levels of thyroid hormone (T3 and T4) may indicate dysfunction of the hypothalamus and pituitary gland. In these case, a high TSH is often produced by a benign

tumour of the pituitary (adenoma). Conversely, low levels of TSH, while blood levels of T3 and T4 are also low, indicates abnormally low function of the pituitary, known as hypopituitarism.

On the other hand, due to the negative feedback described above, abnormally high levels of Thyroid hormone, due to overproduction in the thyroid, results in low TSH levels. This occurs in diseases such as hyperthyroidism or Grave's disease. Conversely, an underproduction of T3 and T4 caused by diseases such as congenital hypothyroidism (cretinism), hypothyroidism or thyroid hormone resistance, gives rise to an increase in the measured TSH.

Clearly both TSH and T3 and T4 should be measured to ascertain where a specific thyroid disfunction is caused by primary pituitary or by a primary thyroid disease.

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
In this method, TSH calibrators, patient specimens and/or controls containing the native antigen are first added to streptavidin coated wells. Biotinylated monoclonal and enzyme labeled antibodies are added and the reactants mixed : these antibodies have high affinity and specificity and are detected against distinct and different epitopes of TSH.

Reaction between the various TSH antibodies and native TSH occurs in the microwells without competition or steric hindrance forming a soluble sandwich complex.

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**AgTSH** = Native Antigen (Variable Quantity)

**EAb(p)** = Enzyme labeled Antibody (Excess Quantity)

**EAb(p)-AgTSH-BtAb(m)** = Antigen-Antibodies Sandwich Complex

**Ka** = Rate Constant of Association

**K-a** = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well though the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

EAb(p)-AgTSH-BtAb(m) + Streptavidina CW. ⇨ Immobilized Complex

**Streptavidin CW.** = Streptavidin immobilized on well  
**Immobilized Complex** = Antibodies-Antigen sandwich bound.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well is quantified by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

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##### 2. Control (1 vial, 1 mL)

Concentration of Control is indicated on the Certificate of Analysis

REF DCE045/1303-0

##### 3. Conjugate (1 vial, 13 mL)

Goat antibody anti TSH conjugated with horseradish peroxidase HRP; mouse antibody anti TSH biotinylated

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Microplates reader (450 nm, 620-630 nm).

#### Note

Store all reagents between 20-8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close immediately after use. Once opened, the microplate is stable until expiry date of the kit.

### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of TSH from 0.2 mIU/L to 20.0 mIU/L.

### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are

needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>6</sub>)

The Calibrators are ready to use, are calibrated against WHO 2<sup>nd</sup> IRP 80/558 and have the following concentrations of TSH:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>	C <sub>6</sub>
mIU/L	0	0.2	0.5	2.5	5.0	10	20

The Calibrators are stable until the expiry date printed on the label. Once opened, the calibrators are stable six months at 2-8°C.

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### 6.3. Preparation of the Sample

TSH assay can be performed in human serum or plasma.

Specimens may be stored at 2-8°C for maximum two days. For longer storage, specimens should be frozen. Frozen specimens should be well mixed after thawing, and before assay. Avoid repeated freezing and thawing.

Serum samples with TSH concentrations greater than 20mIU/L should be diluted 1:2 with C<sub>0</sub>. Correct the result using an appropriate dilution factor.

The Control is ready to use.

### 6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.

- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>6</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C <sub>0</sub> -C <sub>6</sub>	50 µL		
Sample/Control		50 µL	
Conjugate	100 µL	100 µL	
Cover the plate and incubate for 60 minutes at room temperature (22±28°C). To increase the sensibility read the absorbance after 120 min of incubation. Remove the content from each well; wash the wells 6 times with 300 µL of diluted Wash Solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.			
TMB Substrate	100 µL	100 µL	100 µL
Cover the plate and incubate for 20 minutes at room temperature (22±28°C) in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the plate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of TSH for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 8. RESULTS

### 8.1. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C<sub>0</sub>-C<sub>6</sub>) and of each sample.

### 8.2. Calibration curve

Plot the mean value of absorbance (Em) of the calibrators (C<sub>0</sub>-C<sub>6</sub>) against concentration. Draw the best-fit curve through the plotted points. (es:Cubic Spline or Four Parameter Logistic).

### 8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in mIU/L.

## 9. REFERENCE VALUES

Serum samples of apparently healthy women and men were assayed using the TSH ELISA test, with the following results:

	Mean (mIU/L)	Range (mIU/L)
TSH	1.85	0.39 – 6.16

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

## 10. PERFORMANCE AND CHARACTERISTICS

### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (16x) of two different control sera in one assay. The within assay variability is 4.6%.

#### 10.1.2. Inter Assay Variation

Between run variation was determined by replicate measurements (12x) of three different control sera in different lots of kit. The between assay variability is 10.8%.

### 10.2. Accuracy

The recovery of 0.25 – 2.5 – 10 mIU/L of TSH added to sample gave an average value (±SD) of 85.85% ± 2.62% with reference to the original concentrations.

### 10.3. Sensitivity

The lowest detectable concentration of TSH that can be distinguished from the Calibrator 0 is 0.01 mIU/L with a confidence limit of 95%.

### 10.4. Specificity

The cross-reactivity of the thyrotropin ELISA method to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated

by the ratio between the dose of the interfering substance and the amount of TSH needed to obtain the same absorbance.

Thyrotropin (hTSH)	100 %
Thyrotropin - alpha (hTSH- )	9.3 %
Thyrotropin - beta (hTSH- )	152 %
Follicle-Stimulating Hormone (hFSH)	1.0 %
Luteinising hormone (hLH)	1.0 %
Chorionic Gonadotropin (hCG)	<0.1 %

### 10.5. Correlation with CLIA

Diametra TSH ELISA was compared with two others commercially available TSH assays. 54 serum samples were tested.

The two linear regression curves were calculated:

$$(TSH \text{ Diametra}) = 0.96 * (TSH \text{ CLIA}) + 0.06$$

$$r^2 = 0.972$$

$$(TSH \text{ Diametra}) = 0.87 * (TSH \text{ Elisa}) + 0.29$$

$$r^2 = 0.971$$

### 10.6. Hook Effect

Diametra TSH ELISA, a competitive enzyme immunoassay, shows no Hook Effect up to 500 mIU/L.

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Ed. 01/2016

DCM039-10

**DiaMetra S.r.l. Headquater:** Via Calabria 15,  
20090 SEGRATE (MI) Italy

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# FERRITIN ELISA

for routine analysis

Direct immunoenzymatic determination of Ferritin in human serum or plasma

IVD

LOT

See external label

Σ = 96 tests

REF DKO039

## INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of Ferritin concentration in human serum or plasma.

Ferritin ELISA kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

Ferritin is a globular protein found mainly in the liver, which can store about 2'250 iron ( $Fe^{3+}$ ) ions. The ferritin molecule consists of a protein shell (apoferritin) composed of heavy and light subunits, which surrounds a crystalline core containing iron oxide and phosphate.

Ferritin is synthesized in the liver, spleen and numerous other body tissues, with major concentrations found in the liver, spleen, bone marrow, and intestinal mucosa

The ferritin levels measured have a direct correlation with the total amount of iron stored in the body. If ferritin is high there is iron in excess, which would be excreted in the stool. If ferritin is low there is a risk for lack in iron, which sooner or later could lead to anaemia.

In the setting of anaemia, serum ferritin is the most sensitive lab test for iron deficiency anaemia. In contrast, serum ferritin levels are normal or increased in anemia associated with chronic disease. Elevated serum ferritin levels have been observed in acute and chronic liver disease and lymphoid malignancy (leukemia and Hodgkin lymphoma). High serum ferritin levels have also been associated with an elevated risk for myocardial infarction in men. Ferritin is also used as a marker for iron overload disorders, such as haemochromatosis in which the ferritin level may be abnormally raised.

Ferritin is an acute-phase reactant, it is often elevated in the course of disease.

Free iron is toxic to cells, and the body has an elaborate set of protective mechanisms to bind iron in various tissue compartments. Within cells, iron is stored complexed to protein as ferritin or hemosiderin. Apoferritin binds to free ferrous iron and stores it in the ferric state. Under steady state conditions, the serum ferritin level correlates with total body iron stores; thus, the serum ferritin level is the most convenient laboratory test to estimate iron stores.

## 2. PRINCIPLE

Diametra Ferritin ELISA test is based on simultaneous binding of human Ferritin to two monoclonal antibodies, one immobilized on microwell plates and the other conjugated with horseradish peroxidase (HRP).

After incubation the bound/free separation is performed by a simple solid-phase washing.

Then the enzyme HRP in the bound-fraction reacts with the Substrate ( $H_2O_2$ ) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution ( $H_2SO_4$ ) is added.

The colour intensity is proportional to the Ferritin concentration in the sample.

The Ferritin concentration in the sample is calculated based on a standard curve.

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### 3.1. Reagents and materials supplied in the kit

#### 1. Calibrators (6 vials)

CAL0 (3 mL)	REF DCE002/3906-0
CAL1 (1 mL)	REF DCE002/3907-0
CAL2 (1 mL)	REF DCE002/3908-0
CAL3 (1 mL)	REF DCE002/3909-0
CAL4 (1 mL)	REF DCE002/3910-0
CAL5 (1 mL)	REF DCE002/3911-0

#### 2. Control (1 vial, 1 mL)

Control concentration is indicated on the Certificate of Analysis

REF DCE045/3903-0

#### 3. Conjugate (1 vial, 12 mL)

Anti Ferritin antibody conjugated with Horseradish peroxidase (HRP)

REF DCE002/3902-0

#### 4. Coated Microplate (1 breakable microplate)

Anti Ferritin antibody adsorbed on microplate

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Microplates reader (450 nm, 620-630 nm)

#### Note

*Store all reagents between 2-8°C in the dark.*

*Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until the expiry date of the kit.*

*Do not remove the adhesive sheets on the strips unutilised.*

## 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, Calibrators and Controls should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of Ferritin from 5 to 1000 ng/mL.

## 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.

- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are ready to use, are calibrated against the WHO 1<sup>st</sup> IS Ferritin 80/602 and have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
ng/mL	0	5	20	100	400	1000

For sample with concentration over 1000 ng/mL dilute the sample with C<sub>0</sub>.

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### 6.3. Preparation of the Sample

Ferritin determination should be done in human serum or plasma. Specimen can be stored at 2-8°C for at short time (max five days). For longer storage the specimen should be frozen at -20°C. Avoid repeated freezing and thawing. The Control is ready to use.

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Reagent	Calibrator	Sample/Control	Blank
Sample/Control		20 µL	
Calibrator C <sub>0</sub> -C <sub>5</sub>	20 µL		
Conjugate	100 µL	100 µL	
Incubate for 1 hour at room temperature (22-28°C). Remove the content from each well. Wash the wells 3 times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate for 10 minutes at room temperature (22-28°C) in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of Ferritin for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 8. RESULTS

### 8.1. Mean Absorbance

Calculate the mean of the absorbencies (E<sub>m</sub>) corresponding to the single points to the calibration curve (C<sub>0</sub>-C<sub>5</sub>) and of each sample.

### 8.2. Calibration curve

Plot the values of absorbance (E<sub>m</sub>) of the Calibrators (C<sub>0</sub>-C<sub>5</sub>) against concentration. Draw the best-fit curve through the plotted points. (es: Cubic Spline or Four Parameter Logistic).

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## 9. REFERENCE VALUES

The serum values are comprised in the following intervals:

		Mean (ng/mL)	Range (ng/mL)
<b>Women</b>	Premenopausal	53	6 – 180
	Post-menopausal	105	8 – 350
<b>Men</b>		175	20 – 400

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**10.1. Precision***10.1.1. Intra Assay*

Within run variation was determined by replicate (15x) the measurements of three different control sera in one assay. The within assay variability is 7.5%.

*10.1.2. Inter Assay*

Between run variations was determined by replicate (16x) the measurements of three different control sera in different lots. The between assay variability is 6.1%.

**10.2. Specificity**

The cross reaction of the antibody calculated on a weight/weight basis are shown in the table:

Liver Human Iso-Ferritin	100%
Spleen Human Iso-Ferritin	80%
Hearth Human Iso-Ferritin	12%

**10.3. Accuracy**

The recovery of 12.5 – 25 – 50 – 100 – 200 ng/mL of Ferritin added to sample gave an average value ( $\pm$ SD) of 98.66%  $\pm$  2.90%.

The dilution test performed on 3 samples diluted up to 8 times gave an average value ( $\pm$ SD) of 102.11%  $\pm$  5.32%.

**10.4. Sensitivity**

The lowest detectable concentration of Ferritin that can be distinguished from the Calibrator 0 is 0.04 ng/mL at the 95% confidence limit.

**10.5. Correlation with RIA**

Diametra Ferritin ELISA kit was compared to another commercially available Ferritin assay. Serum samples of 22 females and 32 males were analysed according in both test systems.

The linear regression curve was calculated  
 (Ferritin Diametra)= 1.11\*(Ferritin Diasorin) - 10.46  
 $r^2 = 0.972$

**10.6. Hook Effect**

Ferritin ELISA kit, a competitive enzyme immunoassay, shows no Hook Effect up to 50000 ng/mL

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Tel. +39-02-2139184

Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14, 06038 SPELLO (PG)  
Italy

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## CEA ELISA

for routine analysis

Direct immunoenzymatic determination of Carcinoembryonic Antigen (CEA) in human serum or plasma.

IVD

LOT

See external label

Σ = 96 tests

REF DKO051

### INTENDED USE

Diametra CEA ELISA kit is a direct solid phase enzyme immunoassay for quantitative determination of Carcinoembryonic Antigen (CEA) in human serum or plasma.

CEA ELISA kit is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

Carcinoembryonic antigen (CEA) is a glycoprotein, with a molecular weight of 180 kD, involved in cell adhesion. It is normally produced during fetal development, but the production of CEA stops before birth. Therefore, it is not usually present in the blood of healthy adults, although levels are raised in heavy smokers. CEA was identified in human colon cancer tissue extracts. It was later found that serum from individuals with colorectal and other carcinomas had higher levels of CEA than healthy individuals and can be used to monitor the response to colon cancer treatment.

CEA and related genes make up the CEA family belonging to the immunoglobulin superfamily.


The most frequent cancer which causes an increased CEA is cancer of the colon and rectum. Others include cancers of the pancreas, stomach, breast, lung, and certain types of thyroid and ovarian cancer. Benign conditions which can elevate CEA include smoking, infections, inflammatory bowel disease, pancreatitis, cirrhosis of the liver, and some benign tumors in the same organs in which an elevated CEA indicates cancer. Chemotherapy and radiation therapy can cause a temporary rise in CEA due to the death of tumor cells and release of CEA into the blood stream. Benign tumours does not usually cause an increase above 10 ng/mL.

### 2. PRINCIPLE

In this method, the calibrators, the patient specimens and/or the controls (containing the native CEA antigen) are first added to streptavidin coated wells. Biotinylated monoclonal and enzyme labelled antibodies are then added and the reactants mixed: these antibodies have high affinity and specificity and are directed against distinct and different epitopes of CEA.

Reaction between the various CEA antibodies and native CEA occurs in the microwells without competition or steric hindrance, forming a soluble sandwich complex.

The interaction is illustrated by the following equation:

**BtAb(m)** = Biotinylated Monoclonal Antibody (Excess Quantity)

**AgCEA** = Native Antigen (Variable Quantity)

**EAb** = Enzyme labelled Antibody (Excess Quantity)

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**StreptavidinC.W.** = Streptavidin immobilized on well  
**Immobilized complex** = sandwich complex bound to the well

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The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well is quantified by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

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##### 2. Conjugate (1 vial, 13 mL)

Antibodies anti CEA conjugated with Horseradish peroxidase (HRP) and anti CEA biotinylated

**REF** **DCE002/5102-0**

##### 3. Coated Microplate (1 microplate breakable)

Microplate coated with streptavidin

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#### Note

Store all reagents at 20±8°C in the dark.

Open the bag of reagent 3 (Coated Microplate) only when it is at room temperature and close immediately after use; once opened, it is stable until the expiry date of the kit.

### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of CEA from 5.0 to 250 ng/mL.
- CEA has a low clinical sensitivity and specificity as a tumour marker. Clinically an elevated CEA value alone is not of diagnostic value as a test for cancer and should only be used in conjunction with other clinical manifestations (observations) and diagnostic parameters. There are patients with colorectal cancer that do not exhibit elevated CEA values and elevated CEA values do not always change with progression or regression of disease. Smokers demonstrate a higher range of baseline values than non-smokers.

### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems, it is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.

- Maximum precision is required for reconstitution and dispensation of reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are ready to use, are calibrated using a reference preparation which was assayed against the 1<sup>st</sup> IRP 73/601 and have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
ng/mL	0	5	10	25	50	250

Once opened, the Calibrators are stable 6 months at 2-8°C.

### 6.2. Preparation of Wash Solution

Dilute the contents of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

### 6.3. Preparation of the Sample

The specimens shall be human serum or plasma; the usual precautions in the collection of venipuncture samples should be observed.

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For serum preparation, the blood should be collected in a venipuncture tube without additives or anti-coagulants; then allow the blood to clot and centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, they may be stored at temperatures of -20°C for up to 30 days. Avoid repetitive freezing and thawing.

When assayed in duplicate, 0.050 mL of the specimen is required.

For sample with concentration over 250 ng/mL dilute the sample with Calibrator 0.

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Reagent	Calibrator	Sample	Blank
Calibrator C <sub>0</sub> -C <sub>5</sub>	25 µL		
Samples		25 µL	
Conjugate	100 µL	100 µL	
Incubate at room temperature (22±28°C) for 1 hours. Remove the contents from each well. Wash the wells 3 times with 300 µL of diluted Wash Solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22±28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of CEA for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 8. RESULTS

### 8.1. Note

The optical densities (O.D.s) of some calibrators and samples may be higher than 2.0, in such a case, they could be out of the measurement range of the microplate reader. It is therefore necessary, for O.D.s higher than 2.0, to perform a reading at 405 nm (=wavelength of peak shoulder) in addition to 450 nm

(peak wavelength) and 620 (reference filter for the subtraction of interferences due to the plastic).

For microplate readers unable to read the plate at 3 wavelengths at the same time, it is advisable to proceed as follows:

- Read the microplate at 450 nm and at 620 nm.
- Read again the plate at 405 nm and 620 nm.
- Find out the wells whose ODs at 450 nm are higher than 2.0
- Select the corresponding ODs read at 405 nm and multiply these values at 405 nm by the conversion factor 3.0 (where OD 450/OD 405 = 3.0), that is: OD 450 nm = OD 405 nm x 3.0.

Warning: The conversion factor 3.0 is suggested only. For better accuracy, the user is advised to calculate the conversion factor specific for his own reader.

The OD of calibrator 5 should be  $\geq 1.3$ .

### 8.2. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>) and of each sample.

### 8.3. Calibration curve

Plot the mean value of absorbance (Em) of the Calibrators (C<sub>0</sub>-C<sub>5</sub>) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

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## 9. REFERENCE VALUES

Nearly 99% of non-smokers have CEA concentrations less than 5 ng/ml. Similarly 99% of smokers have concentrations less than 10 ng/ml (4).

	ng/mL
<b>Non-smoker</b>	<5 ng/mL
<b>Smoker</b>	<10 ng/mL

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### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (16x) of three different control sera in one assay. The within assay variability is 5.8%.

#### 10.1.2. Inter Assay Variation

Between run variations was determined by replicate measurements (20x) of three different control sera in

different lots of kit. The between assay variability is 9.9%.

### 10.2. Accuracy

The recovery of amounts of 10 - 20 - 40 ng/mL of CEA added to samples gave an average value ( $\pm$ SD) of 100.7%  $\pm$  7.2% with reference to the original concentrations.

### 10.3. Sensitivity

The lowest detectable concentration of CEA that can be distinguished from the Calibrator 0 is 0.11 ng/mL at the 95% confidence limit.

### 10.4. Specificity

In order to assess the specificity of the antibody pair used for the CEA Elisa assay, massive doses of related analytes were spiked in a pool of patient sera:

Analyte	Concentration	Cross Reaction
<b>CEA</b>		100 %
<b>-HCG</b>	250 ng/mL	ND
<b>AFP</b>	5000 ng/mL	ND
<b>CA-125</b>	1000 U/mL	ND
<b>PSA</b>	1000 ng/mL	ND
<b>CA 19-9</b>	1000 U/mL	ND

### 10.5. Correlation

Diametra CEA ELISA was compared to another commercially available CEA assay. 167 serum samples were analysed according in both test systems.

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## CA 15-3 ELISA

for routine analysis

Direct immunoenzymatic determination of CA 15-3 in human serum or plasma

IVD

LOT

See external label

2°C - 8°C

Σ = 96 tests

REF DKO055

### INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of CA 15-3 concentration in human serum or plasma.

CA 15-3 ELISA is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

CA-15-3, is an abbreviation for cancer antigen 15-3. CA 15-3 levels are most useful in following the course of treatment in women diagnosed with breast cancer, especially advanced breast cancer. CA 15-3 levels are rarely elevated in women with early stage breast cancer.

Cancers of the ovary, lung, and prostate may also raise CA 15-3 levels. Elevated levels of CA 15-3 may be associated with non-cancerous conditions, such as benign breast or ovarian disease, endometriosis, pelvic inflammatory disease, and hepatitis. Pregnancy and lactation can also cause CA 15-3 levels to rise.

CA 15-3 is most useful for monitoring patients post-operatively for recurrence, particularly metastatic diseases. 96% of patients with local and systemic recurrence have elevated CA 15-3, which can be used to predict recurrence earlier than radiological and clinical criteria. A 25% increase in the serum CA 15-3 is associated with progression of carcinoma. A 50% decrease in serum CA 15-3 is associated with response to treatment. CA 15-3 is more sensitive than CEA in early detection of breast cancer recurrence. In combination with CA-125, CA 15-3 has been shown to be useful in early detection of relapse of ovarian cancer. CA 15-3 levels are also increased in colon, lung and hepatic tumours.

### 2. PRINCIPLE

CA 15-3 ELISA test is based on simultaneous binding of human CA 15-3 to two monoclonal antibodies, one immobilized on microwell plates, the other conjugated with horseradish peroxidase (HRP).

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Then, the enzyme in the bound-fraction reacts with the Substrate (H<sub>2</sub>O<sub>2</sub>) and the TMB Substrate and develops a blue color that changes into yellow when the Stop Solution (H<sub>2</sub>SO<sub>4</sub>) is added.

The colour intensity is proportional to CA 15-3 concentration in the sample.

CA 15-3 concentration in the sample is calculated using a calibration curve.

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CAL5 REF DCE002/5511-0

2. Control (1 vial, 2 mL)

Concentration of Control is indicated on the Certificate of Analysis REF DCE045/5503-0

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Monoclonal antibody anti CA 15-3 adsorbed on microplate REF DCE002/5503-0

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#### Note

Store all reagents at 2-8°C in the dark.

Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close it

immediately after use; once opened, the microplate is stable until the expiry date of the kit.

#### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservatives. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of CA 15-3 from 0.5 to 240.0 U/mL.

#### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

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#### 6. PROCEDURE

##### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
U/mL	0	15	30	60	120	240

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**CA 15-3 Calibrators and Control have already been pre-diluted and are ready for use.**

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##### 6.3. Preparation of the Sample

The specimens can be human serum or plasma and the usual precautions in the collection of venipuncture samples should be observed.

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Specimen can be stored at 2-8°C for at short time (max five days). For longer storage the specimen should be frozen at -20°C (max 30 days). Avoid repeated freezing and thawing.

**All the samples must be diluted 1:50 with the Serum diluent**, as in the following example:

Serum	20 µL
Serum Diluent	980 µL

Mix gently.

The Control is ready to use.

NB: for samples with concentration over 240 U/mL, dilute further the sample with the Sample diluent (attention: consider this dilution in the calculation of final result).

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Reagent	Calibrator	Samples/ Control	Blank
Control / Diluted samples		200 µL	
Calibrator C <sub>0</sub> -C <sub>5</sub>	200 µL		
Incubate for 1 hour at 37°C without covering the plate. Remove the content from each well. Wash the wells 5 times with 300 µL of diluted Wash Solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.			
Conjugate	200 µL	200 µL	
Incubate for 1 hour at 37°C without covering the plate. Remove the content from each well. Wash the wells 5 times with 300 µL of diluted Wash Solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22±28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake gently the microplate. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

#### 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of CA 15-3 for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

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##### 8.2. Calibration curve

Plot the mean value of absorbance (E<sub>m</sub>) of the Calibrators (C<sub>0</sub>-C<sub>5</sub>) against concentration. Draw the best-fit curve through the plotted points. (Es: Four Parameter Logistic).

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### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (15x) of three different control sera in one assay. The within assay variability is 7.5%.

#### 10.1.2. Inter Assay Variation

Between run variations was determined by replicate measurements (10x) of five different control sera in different lots of kit. The between assay variability is 11.4%.

### 10.2. Accuracy

The recovery on three serum samples spiked with 25 - 50 - 100 U/mL of antigen gave an average value ( $\pm$  SD) of 106.61%  $\pm$  8.86%.

The dilution test performed on three sera diluted 2 - 4 - 8 - 16 times gave an average value ( $\pm$ SD) of 103.71%  $\pm$  8.56%.

### 10.3. Sensitivity

The lowest detectable concentration of CA 15-3 that can be distinguished from the Calibrator zero is 0.40 U/mL at the 95% confidence limit.

### 10.4. Specificity

The cross reaction of the antibody are shown in the table:

Antigens	Concentration	% Cross-reactivity
CA 15-3	---	100.0%
CA 125	1,000 U/mL	N.D.
CA 19-9	1,000 U/mL	N.D.
PSA	1,000 ng/mL	N.D.
PAP	1,000 ng/mL	N.D.
AFP	10,000 ng/mL	N.D.
CEA	1,000 ng/mL	N.D.

### 10.5. Correlation

Diametra CA 15-3 ELISA kit was compared with a commercial reference method. 54 serum samples were analysed according in both test system. The linear regression curve was calculated:

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## CA19-9 ELISA

for routine analysis

Direct immunoenzymatic determination of CA19-9 in human serum or plasma.

IVD

LOT

See external label

Σ = 96 tests

REF DKO056

### INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of CA19-9 concentration in human serum or plasma.

CA19-9 ELISA is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

CA19-9 represents the most important tumour marker of the mucin type glycoproteins group, that are associated to gastrointestinal cancer pathology.

The immunohistologic distribution of CA19-9 in tissues is consistent with the quantitative determination of higher CA19-9 concentrations in cancer than in normal or inflamed tissues. Recent reports indicate that the serum CA19-9 level is frequently elevated in the serum of subjects with various gastrointestinal malignancies, such as pancreatic, colorectal, gastric and hepatic carcinomas. Together with CEA, elevated CA19-9 is suggestive of gall bladder neoplasm in the setting of inflammatory gall bladder disease. This tumour-associated antigen may also be elevated in some non-malignant conditions.

It has been shown that a persistent elevation in serum CA19-9 value following treatment may be indicative of occult metastatic and/or residual disease. Research studies demonstrate that serum CA19-9 values may have utility in monitoring subjects with the above-mentioned diagnosed malignancies.

A persistently rising serum CA19-9 value may be associated with progressive malignant disease and poor therapeutic response. A declining CA19-9 value may be indicative of a favourable prognosis and good response to treatment.

### 2. PRINCIPLE

CA19-9 ELISA test is based on simultaneous binding of human CA19-9 to two monoclonal antibodies, one immobilized on microwell plates, the other conjugated with horseradish peroxidase (HRP). After incubation, the separation bound-free is obtained with a simple solid-phase washing. Then, the enzyme in the bound-fraction reacts with the Substrate (H<sub>2</sub>O<sub>2</sub>) and the TMB Substrate and develops a blu color that changes into yellow when the Stop Solution (H<sub>2</sub>SO<sub>4</sub>) is added.

The colour intensity is proportional to CA19-9 concentration in the sample. CA19-9 concentration in the sample is calculated using a calibration curve.

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##### 2. Control (1 vial, 1 mL)

The concentration of the Control is lot-specific and is indicated on the Certificate of Analysis

REF DCE045/5603-0

##### 3. Conjugate (1 vial, 12 mL)

Antibody anti CA19-9 conjugated with horseradish peroxidase (HRP)

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Antibody anti CA19-9 adsorbed on microplate

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#### Note

Store all reagents at 20±8°C in the dark.

Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until the expiry date of the kit. Do not remove the adhesive sheets on the strips unutilised

#### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents contain small amounts of Proclin 300<sup>R</sup> as preservative. Avoid the contact with skin or mucosa.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the determination of CA19-9 from 0.2 to 240.0 U/mL.

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#### 6. PROCEDURE

##### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are ready to use and have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
U/mL	0	15	30	60	120	240

Once opened, the Calibrators are stable 6 months at 2-8°C

##### 6.2. Preparation of Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C. In concentrated wash solution is possible to observe the presence of crystals, in this case mix at room temperature until complete dissolution of the crystals; for greater accuracy dilute the whole bottle of concentrated wash solution to 500 mL taking care also to transfer the crystals, then mix until crystals are completely dissolved.

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Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, the samples may be stored at -20°C for up to 30 days. Avoid repetitive freezing and thawing.

For sample with concentration over 240 U/mL, dilute the sample with the Calibrator 0.

The Control is ready to use.

##### 6.4. PROCEDURE

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test

results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Sample/ Control		100 µL	
Calibrator C <sub>0</sub> -C <sub>5</sub>	100 µL		
Incubation 1 h at +37°C. Remove the contents from each well; wash the wells 3 times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
Conjugate	100 µL	100 µL	
Incubation 1 h at +37°C. Remove the contents from each well; wash the wells 3 times with 300 µL of diluted wash solution. <b>Washing:</b> follow the same indications of the previous point.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature 22÷28°C for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of CA19-9 for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

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### 8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in U/mL.

## 9. REFERENCE VALUES

	CA19-9
Healthy people (male and non-pregnant females)	35 U/mL

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### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate (10x) the measurement of three different control sera in one assay. The within assay variability is ≤ 3.4%.

#### 10.1.2. Inter Assay Variation

Between run variation was determined by replicate (10x) the measurement of three different control sera with different lots of kit. The between assay variability is ≤ 6.6%.

### 10.2. Sensitivity

The lowest detectable concentration of CA19-9 that can be distinguished from the Calibrator zero is 0.18 U/mL at the 95 % confidence limit.

### 10.3. Specificity

The cross-reaction of the CA19-9 assay was assessed by measuring apparent dose of various potentially cross-reagents. The results are shown in the following table:

Antigens	Concentration	% Cross-reactivity
CA19-9	---	100.0%
CA 15-3	1,000 U/mL	3.0%
CA 125	400 U/mL	Not Detected
PSA	100 ng/mL	Not Detected
PAP	1,000 ng/mL	Not Detected
AFP	10,000 ng/mL	Not Detected
CEA	250 ng/mL	Not Detected

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#### 10.4. Dilution

Three patient serum samples containing concentrations of CA19-9 up to 100 IU/mL, were serially diluted with Calibrator zero provided in the kit up to a dilution of 8 and assayed.

The percentage of the mean of recovered dose are 112% (s.d. 3.1%).

#### 10.5. Recovery

Three normal serum samples were spiked with 200, 100, 50, 25 and 0 U/mL of CA19-9 and then assayed.

The percentage of the mean of recovered dose is 92.4% (s.d. 8.8%).

#### 10.6. Correlation

CA19-9 ELISA assay was compared to an established method using Roche Modular; 22 serum samples were assayed with both methods.

A linear regression analysis was performed and the following results were obtained:

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## Anti TG ELISA

for routine analysis


Quantitative determination of antibodies against Thyroglobulin (TG) in human serum or plasma

IVD

LOT

See external label

Σ = 96 tests

REF DKO115

### INTENDED USE

Indirect Immunoenzymatic colorimetric method for quantitative determination of Anti-Tg concentration in human serum or plasma.

Anti TG ELISA kit is intended for laboratory use only.

### 1. CLINICAL SIGNIFICANCE

Thyroglobulin (Tg) is a protein produced by and used entirely within the thyroid gland. Thyroglobulin is used by the thyroid gland to produce the thyroid hormones thyroxine (T4) and triiodothyronine (T3).

Thyroglobulin synthesis is stimulated at the transcriptional level by thyroid-stimulating hormone (TSH). Thyroglobulin (Tg) is a well-known target for autoantibodies occurring in thyroid autoimmunity (Graves' disease and Hashimoto's thyroiditis). Anti-Tg antibodies mostly belong to the IgG class. Low to moderate levels of anti-Tg antibodies can be found in sera of other autoimmune patients (eg systemic lupus erythematosus or Sjogren syndrome). In some cases anti-Tg positive sera may show negativity for other type of anti-thyroid antibodies - anti-TPO.

Therefore, combined determination of both types of anti-thyroid antibodies (anti-TPO + anti-Tg) provides most sensitive laboratory diagnostic tool for thyroid autoimmunity. Separately from autoimmunity, anti-Tg antibodies may develop in patients suffering from thyroid cancer. High level of anti-Tg in such patients may interfere with correct determination of serum thyroglobulin, which serves as tumour marker for therapy control in this group of patients.

### 2. PRINCIPLE

This test is based on a two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated with the antigen. Antibodies from the specimen bind to the coated antigen on the microwell surface. Unbound material is removed by the washing procedure. A second antibody directed towards human IgG and labeled with peroxidase enzyme is then added into the microwells. After the subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of a chromogen substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of specific antibodies in the specimen.

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- Calibrators** (6 vials, 1.2 mL each)  
Phosphate buffer 0.1M, NaN<sub>3</sub> <0.1%, human serum  
CAL0 **REF DCE002/11506-0**  
CAL1 **REF DCE002/11507-0**  
CAL2 **REF DCE002/11508-0**  
CAL3 **REF DCE002/11509-0**  
CAL4 **REF DCE002/11510-0**  
CAL5 **REF DCE002/11511-0**
- Controls** (2 vials, 1.2 mL each, ready to use)  
Phosphate buffer 0.1M, NaN<sub>3</sub> <0.1%, human serum  
Negative Control **REF DCE045/11501-0**  
Positive Control **REF DCE045/11502-0**
- Sample Diluent** (1 vial, 100 mL)  
Phosphate buffer 0.1M, NaN<sub>3</sub> <0.1% **REF DCE053-0**
- Conjugate** (1 vial, 15 mL)  
Anti h-IgG conjugated with horseradish peroxidase (HRP), BSA 0.1%, Proclin <0.0015% **REF DCE002/11502-0**
- Coated Microplate** (1 breakable microplate)  
Microplate coated with antigenic Thyroglobulin **REF DCE002/11503-0**
- TMB Substrate** (1 vial, 15 mL)  
H<sub>2</sub>O<sub>2</sub>-TMB (0,26 g/L) (avoid any skin contact) **REF DCE004/-0**
- Stop Solution** (1 vial, 15 mL)  
0.15M Sulfuric acid (avoid any skin contact) **REF DCE005-0**
- 50X Conc. Wash Solution** (1 vial, 20 mL)  
NaCl 45 g/L; Tween-20 55 g/L **REF DCE006-0**

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Microplates reader (450 nm, 620-630 nm)

#### Notes

Store all reagents between 20±8°C in the dark.

Open the bag of reagent 5 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until the expiry date of the kit.

#### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
- Some reagents contain small amounts of Sodium Azide (NaN<sub>3</sub>) or Proclin 300<sup>R</sup> as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.

#### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips.**

For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.

- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

#### 6. PROCEDURE

##### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The assay system is calibrated in relative arbitrary units. The Calibrators are ready to use and have the following concentration:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
AU/mL	0	2	4	8	32	128

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##### 6.2. Preparation of the Sample

For determination of Anti-Tg human serum or plasma are the preferred sample matrix.

**Serum samples have to be prediluted with sample diluent 1:100;** for example, 10 µL of sample may be diluted with 990 µL of sample diluent.

The patients need not to be fasting, and no special preparations are necessary. Collect blood by venipuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Samples may be stored refrigerated at 2-8°C for at least 5 days. For longer storage of up to six months samples should be stored frozen at -20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

Do not use microbiologically contaminated samples, as highly lipemic or hemolyzed samples.

The Controls are ready to use.

##### 6.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (50X) with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

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Reagent	Calibrator	Sample or Controls	Blank
Calibrator C <sub>0</sub> -C <sub>5</sub>	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
Conjugate	100 µL	100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution. <b>Washing:</b> follow the same indications of the previous point.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. RESULTS

### 7.1. Calibration curve

For Anti TG ELISA kit a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable. We recommend using a Lin-Log curve.

First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may also be connected with straight line segments. The concentration of unknowns may then be estimated from the calibration curve by interpolation.

### Typical Results (example only)

The figures below show typical results for Anti-Tg. These data are intended for illustration only and should not be used to calculate results from another run.

CAL	OD1	OD2	mean	C1	C2	mean	CV%
0	0.025	0.032	0.0285	0.045	0	0.0225	141,42
1	0.231	0.24	0.2355	1.88	2.12	2	8,49
2	0.429	0.442	0.4355	3.96	4.2	4.08	4,16
3	0.767	0.759	0.763	8.2	7.77	7.985	3,81
4	1.884	2.065	1.9745	31.9	32.36	32.13	1,01
5	3.190	3.316	3.253	128.2	127.5	127.85	0,39

## 8. REFERENCE VALUES

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti-Tg tests:

	Anti TG [AU/mL]
Normal	< 4
Elevated	4

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Positive results should be verified concerning the entire clinical status of the patient. Also every decision for therapy should be taken individually.

It is recommended that each laboratory establishes its own normal and pathological ranges of serum Anti-Tg.

## 9. PERFORMANCE AND CHARACTERISTICS

### 9.1. Specificity

Comparison test against a commercial reference kit, performed on 66 sera (27 of them positive sera and 39 negative sera) shows a 94.9% specificity.

### 9.2. Sensibility

Comparison test against a commercial reference kit, performed on 66 sera (27 of them positive sera and 39 negative sera) shows a 96.3% sensibility.

### 9.3. Detection limit

The lowest concentration of anti Tg that can be distinguished from Calibrator 0 is about 0.11 AU/mL with a confidence limit of 95%.

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**Ed. 01/2015**

**DCM115-5**

**DiaMetra S.r.l. Headquarter:** Via Calabria 15,  
20090 SEGRATE (MI) Italy  
Tel. +39-02-2139184  
Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14  
06038 SPELLLO (PG) Italy  
Tel. +39-0742-24851  
Fax +39-0742-316197  
E-mail: [info@diametra.com](mailto:info@diametra.com)

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# Anti TPO ELISA

for routine analysis

Quantitative determination of antibodies against Thyroperoxidase in human serum or plasma

IVD

LOT

See external label

$\Sigma = 96$  tests

REF DKO116

## INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of Anti-TPO concentration in human serum or plasma.

Anti TPO ELISA kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

Anti TPO antibodies (formerly thyroid microsomal antibodies) are directed against a target protein - thyroid peroxidase (TPO) - located in the smooth endoplasmic reticulum of thyroid cells.

Thyroid peroxidase or Thyroperoxidase (TPO) is an enzyme mainly expressed in the thyroid that liberates iodine for addition onto tyrosine residues on thyroglobulin for the production of thyroxine (T4) or triiodothyronine (T3) (thyroid hormones).

It is a frequent epitope of autoantibodies in autoimmune thyroid disease. In autoimmune disorders, your immune system makes antibodies that mistakenly attack normal tissue such as thyroid peroxidase. Antibodies that attack the thyroid gland cause inflammation and impaired function of the thyroid.

The presence of anti-TPO antibodies in serum is associated with thyroid autoimmune diseases (Graves' disease and Hashimoto's thyroiditis). Anti-TPO antibodies mostly belong to the IgG class.

Low to moderate levels of serum anti-TPO antibodies can be found in some other autoimmune pathology (eg systemic lupus erythematosus or Sjogren syndrome) and, rarely, in apparently healthy subjects (especially elderly women). Anti-TPO antibodies are more sensitive in diagnostics of thyroid autoimmune diseases than anti-thyroglobulin (anti-TG) antibodies.

However, in some cases anti-TG positive sera may be negative for anti-TPO. Therefore, combined determination of both types of anti-thyroid antibodies (anti-TPO + anti-TG) provides a more sensitive laboratory diagnostic tool for thyroid autoimmunity.

## 2. PRINCIPLE

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by the antigen. Antibodies from the specimen bind coated antigen on the microwell surface. Unbound material is removed by washing procedure. Second antibodies directed towards human Ig and labeled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogensubstrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of specific antibodies in the specimen.

## 3. REAGENTS, MATERIALS AND INSTRUMENTATION

### 3.1. Reagent and material supplied in the kit

- Calibrators** (6 vials, 1.2 mL each)  
Phosphate buffer 0.1M,  $\text{NaN}_3 < 0.1\%$ , human serum  
CAL0 **REF DCE002/11606-0**  
CAL1 **REF DCE002/11607-0**  
CAL2 **REF DCE002/11608-0**  
CAL3 **REF DCE002/11609-0**  
CAL4 **REF DCE002/11610-0**  
CAL5 **REF DCE002/11611-0**
- Controls** (2 vials, 1.2 mL each, ready to use)  
Phosphate buffer 0.1M,  $\text{NaN}_3 < 0.1\%$ , human serum  
Negative Control **REF DCE045/11601-0**  
Positive Control **REF DCE045/11602-0**
- Sample Diluent** (1 vial, 100 mL)  
Phosphate buffer 0.1M,  $\text{NaN}_3 < 0.1\%$   
**REF DCE053-0**
- Conjugate** (1 vial, 15 mL)  
Anti h-IgG conjugated with peroxidase, BSA 0.1%, Proclin  $< 0.0015\%$   
**REF DCE002/11602-0**
- Coated Microplate** (1 breakable microplate)  
Microplate coated with antigenic Thyroperoxidase)  
**REF DCE002/11603-0**
- TMB Substrate** (1 vial, 15 mL)  
 $\text{H}_2\text{O}_2$ -TMB (0,26 g/L) (avoid any skin contact)  
**REF DCE004-0**
- Stop Solution** (1 vial, 15 mL)  
Sulphuric acid 0.15M (avoid any skin contact)  
**REF DCE005-0**
- 50X Conc. Wash Solution** (1 vial, 20 mL)  
NaCl 45 g/L; Tween-20 55 g/L **REF DCE006-0**

### 3.2. Reagents necessary not supplied

Distilled water.

### 3.3. Auxiliary materials and instrumentation

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#### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of Calibrators and controls for this product has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrator and the Controls should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of Sodium Azide (NaN<sub>3</sub>) or Proclin 300<sup>R</sup> as preservatives. Avoid the contact with skin or mucosa.
- Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, allow scroll through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.

#### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding

contamination of the conjugate; this procedure is highly recommended when the kit is processed using analyzers which are not equipped with disposable tips. For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.

- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated should not be used in the assay. Highly lipemic or haemolysed specimens should similarly not be used
- Plate readers measure vertically. Do not touch the bottom of the wells.

#### 6. PROCEDURE

##### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The assay system is calibrated in relative arbitrary units. The Calibrators have the following concentration:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
AU/mL	0	5	10	20	80	320

##### 6.2. Preparation of the Sample

For the determination of anti TPO antibodies human serum or plasma is the preferred sample matrix.

**Serum samples have to be prediluted with sample diluent 1:100;** for example 10 µL of sample may be diluted with 990 µL of sample diluent.

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Samples may be stored refrigerated at 2 - 8°C for at least 5 days. For longer storage of up to six months samples should be stored frozen at -20°C. To avoid repeated thawing and freezing the samples should be aliquoted.

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The Controls are ready to use.

##### 6.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (50x) with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

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Reagent	Calibrator	Sample or Controls	Blank
Calibrator C <sub>0</sub> -C <sub>5</sub>	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the contents from each well, wash the wells three times with 300 µL of diluted wash solution. <b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. <b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times.			
Conjugate	100 µL	100 µL	
Incubate 30 minutes at room temperature (22-28°C). Remove the contents from each well, wash the wells three times with 300 µL of diluted wash solution. <b>Washing:</b> follow the same indications of the previous point.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22-28°C).			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

## 7. RESULTS

### 7.1. Calibration curve

For Anti-TPO a 4-Parameter-Fit with lin-log coordinates for optical density and concentration is the data reduction method of choice. Smoothed-Spline Approximation and log-log coordinates are also suitable. However we recommend using a Lin-Log curve.

First calculate the averaged optical densities for each calibrator well. Use lin-log graph paper and plot the averaged optical density of each calibrator versus the concentration. Draw the best fitting curve approximating the path of all calibrator points. The calibrator points may also be connected with straight line segments. The

concentration of unknowns may then be estimated from the calibration curve by interpolation.

### Typical results (example only)

The figures below show typical results for Anti-TPO. These data are intended for illustration only and should not be used to calculate results from another run.

CAL	OD1	OD2	mean	C1	C2	Mean	CV%
0	0,012	0,013	0,013	0	0,03	0,02	141,42
1	0,125	0,121	0,123	5,65	5,45	5,55	2,49
2	0,212	0,210	0,211	9,94	9,84	9,89	0,71
3	0,395	0,422	0,409	19,37	20,83	20,10	5,12
4	1,303	1,221	1,262	84,12	76,43	80,28	6,77
5	2,521	2,474	2,498	330,2	310,9	320,52	4,26

## 8. REFERENCE VALUES

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti-TPO tests:

	Anti TPO [AU/mL]
Normal	< 20
Elevated	20

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It is recommended that each laboratory establishes its own normal and pathological ranges of serum Anti-TPO.

## 9. PERFORMANCE AND CHARACTERISTICS

### 9.1. Specificity

Comparison test against a commercial reference kit, performed on 78 sera (50 of them positive sera and 28 negative sera) shows a 92.9% specificity.

### 9.2. Sensitivity

Comparison test against a commercial reference kit, performed on 78 sera (50 of them positive sera and 28 negative sera) shows a 94.0% sensibility.

### 9.3. Detection limit

The lowest concentration of anti-TPO that can be distinguished from zero Calibrator is 0.09 AU/mL with a confidence limit of 95%.

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Ed. 01/2015

DCM137-2

**DiaMetra S.r.l. Headquater:** Via Calabria 15,  
20090 SEGRATE (MI) Italy  
Tel. +39-02-2139184  
Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14  
06038 SPELLLO (PG) Italy  
Tel. +39-0742-24851  
Fax +39-0742-316197  
E-mail: [info@diametra.com](mailto:info@diametra.com)

DCM137-2  
Ed. 01/2015

## Total PSA

for routine analysis

Enzyme immunoassay for the quantitative determination of total PSA (Prostate Specific Antigen) in human serum and plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 test

REF DKO137

### INTENDED USE

DiaMetra Total PSA ELISA is used for the quantitative determination of total Prostate Specific Antigen (t-PSA) in human serum or plasma samples. The determination of total PSA levels is used to estimate the risk of prostate carcinoma in men in conjunction with digital rectal examination (DRE) or to monitor the effectiveness of prostate carcinoma treatment in patients.

### 1. CLINICAL SIGNIFICANCE

Prostate cancer is the most frequent type of cancer found in man and is the second cause of death due to cancer in males. Until recently, digital rectal examination (DRE) was frequently used as only diagnostic modality for the detection of early stages of prostate cancer. In the recent years the determination of serum PSA levels has become the most accepted method to improve the diagnostic specificity of DRE. Although PSA is a tissue specific protein and is not solely tumor specific, it has become the most important marker for prostate carcinoma, showing a better specificity than other biochemical markers used in this context (PAP, total alkaline phosphatase, carcinoembryonic antigen, etc.).

In 1979, Wang et al isolated a specific antigen for normal prostate tissue and called this protein PSA. As demonstrated by immunohistological studies, PSA is localized in the cytoplasm of prostate acinar cells, ductal epithelium and in the secretion on the ductal lumina, present in normal, benign hyperplastic and malignant prostate tissues as well metastatic prostate cancer and in seminal plasma. If the structural integrity of the prostate is disturbed and/or the gland size is increased, the amount of PSA in the blood plasma may become elevated. An elevation of PSA levels to values higher than 3-4 ng/ml has been reported for patients with either benign prostatic hypertrophy (BPH) or prostate carcinoma. At this threshold follow-up examinations that allow to differentiate between these two conditions are recommended.

The determination of PSA serum levels is not only important for the screening of patients for prostate cancer, but also for monitoring patients that have been treated for this disease. Here regular PSA measurements are an important tool to examine the potential and actual effectiveness of surgery or other therapies. An increase of PSA in patients after radical prostatectomy or radiotherapy may allow an earlier discovery of residual or recurrent carcinoma.

### 2. PRINCIPLE OF THE METHOD

This assay is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the sandwich principle.

The microtiter wells are coated with an antibody, directed towards an epitope of an antigenic PSA. An aliquot of patient serum is incubated in the coated well with enzyme conjugated second antibody (E-Ab), directed towards a different region of the antigen molecule. After incubation the unbound E-Ab is washed off. The amount of bound E-Ab is proportional to the concentration of antigen in the sample. After adding the substrate solution, the intensity of colour developed is proportional to the antigen concentration in the sample. The measured ODs of the Calibrators are used to construct a calibration curve against which the unknown samples are calculated.

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Contain non-mercury preservative.

##### 2. Controls (2 vials, 0.5 mL each)

For concentration see the Certificate of Analysis.

Contains non-mercury preservative.

Low Control REF DCE045/13701-0  
High Control REF DCE045/13702-0

##### 3. Coated Microplate (1 breakable microplate)

Monoclonal antibody anti PSA adsorbed on microplate

REF DCE002/13703-0

##### 4. Enzyme Conjugate (1 vial, 12 mL)

Monoclonal antibody anti PSA conjugated to horseradish peroxidase. Contains non-mercury preservative.

REF DCE002/13702-0

##### 5. Substrate Solution (1 vial, 12 mL)

Contains TMB (tetramethylbenzidine) and H<sub>2</sub>O<sub>2</sub>

REF DCE004/13704-0

##### 6. Stop Solution (1 vial, 14 mL)

Contains sulphuric acid (H<sub>2</sub>SO<sub>4</sub>) 0.5M (avoid contact with the Stop Solution, it may cause skin irritations and burns)

REF DCE005/13705-0

### 3.2. Material required but not provided

- Distilled water.
- Precision micropipettes (volume: 25 µL and 100 µL) with disposable tips.
- ELISA photometer with 450 nm and 630 nm filters.
- Timer with 60 min. range or higher.
- Microplate washer (optional).
- Vortex or similar mixing tools.
- Container for the proper handling of waste and samples after use.

### 4. STORAGE AND STABILITY

- Store the kit and components at +2 to + 8°C
- Bring to room temperature (18-25°C) at least 30 minutes before use. After use put back into the refrigerator. Avoid long time storage at room temperature.
- Do not use the kit or components after the expiry date. For expiry date of the original packed kit see kit label.
- Close the bottles immediately after use.
- Store the plate incl. desiccant in the provided zip-lock pouch. Modules that are not used should always be stored under this condition.
- Ensure that kit components do not freeze.

### 5. PRECAUTIONS

- ELISA kits are only for in vitro diagnostic use by professionals.
- Serum and plasma samples should be treated as potentially infectious materials. Wear gloves and proper laboratory attire when handling sample materials. Do not eat, drink or smoke in areas where specimen or kit reagents are handled. Do not pipette with the mouth. In case of skin contact, wash with a germicidal soap and copious amounts of water. Seek medical advice if indicated.
- The PSA Calibrators and controls are of human origin. They have been tested and confirmed negative for HIV, HBsAg and HCV. However, all Calibrators should be treated as potential biohazards.
- Due to the potentially infectious character of samples and kit components all materials that have come in contact with these materials should be sterilized and disposed of according to local regulations. This also includes the liquid waste.
- The assay reagents contain preservatives, TMB, H<sub>2</sub>O<sub>2</sub> or sulphuric acid and may be harmful if ingested. A direct skin or mucosa contact should be avoided. In case of skin contact, wash thoroughly with water and seek medical attention if required.
- The stop solution contains H<sub>2</sub>SO<sub>4</sub>. Since the H<sub>2</sub>SO<sub>4</sub> used to terminate the color reaction is corrosive, the instrumentation employed to dispense it should be thoroughly cleaned after use.
- Do not interchange reagents from different lots or different suppliers.
- Avoid reagent or sample carry-over by using fresh tips for solutions and samples.
- Do not use test kit if zip lock pouch or bottles have been damaged.

## 6. GUIDELINE FOR SAMPLE COLLECTION, PREPARATION AND STORAGE

### 6.1. Sample collection

Blood samples are collected by venipuncture. As different factors could influence the PSA level in blood, doctors should ensure that the patient has avoided the following conditions before taking the blood sample.

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- biking
- sexual intercourse (ejaculation)
- Manipulation of the prostate during medical examinations like DRE, transrectal prostatic ultrasound etc.
- Prostatitis
- Liver dysfunction

The following conditions may lead to a decrease of PSA levels:

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### 6.3. Storage of samples

For the assay either fresh serum or plasma samples can be used.

If not used immediately they can be stored at 2-8°C for 1 week. In case of longer storage, freeze at -20°C. A repeated freezing and thawing of samples should be avoided.

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## 7. ASSAY PROCEDURE

### 7.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are calibrated against the WHO 96/670 and have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
ng/mL	0	1.56	3.12	6.25	12.50	25.00

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- Prior to use bring all reagents, Calibrators, controls, and samples to room temperature (18-25°C).
- Check that all components are not expired and take care that bottles and plate (inclusive pouch) are not damaged.
- Format the required microplate wells. Keep in mind that all measurements should per performed as duplicate. Document position of wells and respective samples, Calibrators and controls to ensure later identification.

Put any unused microwell modules back into the zip lock bag with the desiccant, seal bag and store at 2-8 °C.

- Samples with an expected PSA value higher than 25 ng/mL should be diluted with the sample diluent.

Reagent	Calibrator	Sample	Control
Calibrator C0-C5	25 µL		
Control			25 µL
Sample		25 µL	
Incubate 5 minutes at room temperature (18-25°C).			
Enzyme Conjugate	100 µL	100 µL	100 µL
Mix by moving plate on the table (10 sec) Incubate 1 hour at room temperature (18-25°C). Remove solution from the wells by aspirating the liquid or by decanting it. If decanting, tap plate on adsorbent paper to remove residual liquid. For washing fill plate with distilled water and wait 15 sec before removing the distilled water; wash 5x to 6x. <u>We recommend the following procedure:</u> wash wells 6-times with 250 µL / well distilled water. Preferably use an automated washing procedure. If washing manually take care that the washing solution remains in each well for the same time. This is necessary to receive lowest possible CV-values.			
Substrate Solution	100 µL	100 µL	100 µL
Incubate 20 minutes at room temperature (18-25°C) in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read absorbencies (OD) at 450 nm (blanking 630 nm)			

## 8. CALCULATION OF RESULTS

1. Calculate the average absorbance values for each set of Calibrators, controls and patient samples.
2. Using linear graph paper, construct a calibration curve by plotting the mean absorbance obtained from each Calibrator against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the calibration curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred method. Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this calibration curve. Samples with concentrations higher than that of the highest Calibrator have to be further diluted or reported as > 25 ng/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

### 8.1. Example of Typical Calibration Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Calibrator	Optical Units (450 nm)
Calibrator 0 (0 ng/mL)	0.05
Calibrator 1 (1,56 ng/mL)	0.24
Calibrator 2 (3,12 ng/mL)	0.39
Calibrator 3 (6,25 ng/mL)	0.74
Calibrator 4 (12.5 ng/mL)	1.27
Calibrator 5 (25,0 ng/mL)	2.01

## 9. QUALITY CONTROL

- It is recommended that internal controls are used in every assay in duplicate. Control results should be within established ranges and should preferably represent low, medium, and high concentrations.
- The risk for the patient is mainly due to false negative results around the reschild (cut-off) value of PSA < 4.0 ng/mL. It is therefore highly recommended to validate the kit and laboratory via external trials (e.g. DGKC).

## 10. EXPECTED VALUES AND LIMITATIONS OF THE PROCEDURE

The generally recommended threshold for follow-up examinations is:

Cut-off value: 3.0-4.0 ng PSA /mL

Healthy men generally have a PSA concentration lower than 4.0 ng/mL. If the PSA concentration is equal or higher than 4.0 ng/ml follow-up examinations are highly recommended. This PSA concentration indicates an elevated risk for prostate cancer but might also be caused by BPH.

Please note that that the 4 ng/ml threshold is only a guideline value. In the literature it is discussed that modifications according to age and ethnological background might be useful e.g. that for younger men the threshold should be lower than for older men. If possible, it is recommended for each laboratory to establish its own specific values that take into consideration a population indigenous to the area where the laboratory is located.

It is important to keep in mind that some prostate tumors do not cause elevated PSA levels so that PSA measurements should never replace DRE but should only be used in conjunction with DRE.

As elevated PSA levels might also be caused by non-cancerous conditions follow examinations might try to increase the diagnostic specificity of t-PSA values. In the literature PSA density, PSA velocity and the ratio of f-PSA to t-PSA are discussed to improve discrimination between cancerous and non-cancerous conditions and might be used to reduce unnecessary prostate biopsies. But only a prostate biopsy can finally show if a prostate carcinoma is present or not.

**Note: PSA values can only be used to estimate the cancer risk. They should always be interpreted in conjunction with other clinical findings and should not be used as a sole basis for prostate cancer diagnosis.**

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Detection limit

The limit of detection for this kit is 0.2 ng/mL.

### 11.2. Precision

Intra- and inter-assay precisions were established by analysing three patient sera of different PSA concentrations. The results are shown in Tables 1 and 2.

Table 1. Intra-assay precision

Patients	Number of replicates	Mean ng/mL	SD ng/mL	CV %
1	24	12.52	0.65	6.0
2	24	3.44	0.13	3.9
3	32	0.83	0.07	8.8

Table 2. Inter-assay precision

Patients	Number of replicates	Mean ng/mL	SD ng/mL	CV %
1	4	12.28	0.82	6.7
2	4	3.33	0.266	7.98

### 11.3. Recovery

A known amount of PSA was added to three patient sera and the quantities recovered were measured. The results are shown in Table 3.

Table 3. Recovery

Sample	Expected value (ng/mL)	Observed value (ng/mL)	Recovery %
1	6.30	6.40	102
2	4.67	4.56	98
3	10.10	10.91	108

### 11.4. Specificity

The antibodies used in this kit are highly specific for total PSA (free PSA & PSA-ACT-complex), with a relatively low cross-reactivity to other proteins and polypeptides, lipids or chemotherapeutic agents that might be present in patient samples.

Table 4. Specificity

Antigens	Amount added	Cross reaction
<b>Proteins</b>		
AFP	10 µg/mL	No
CEA	10 µg/mL	No
HCG	10 µg/mL	No
Lactalbumin	10 µg/mL	No
PAP	1 µg/mL	No
<b>Interfering substances</b>		
Bilirubin	0.2 mg/mL	No
Hemoglobin*	0.1 mg/ml	No
Triglyceride	15 mg/mL	No
<b>Chemotherapeutic Agents</b>		
Cyclophosphamide	800 µg/mL	No
Doxorubicin * HCl	20 µg/mL	No
Diethylstilbestrol	2 µg/mL	No
Flutamide	10 µg/mL	No
Methotrexate	50 µg/mL	No

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### 11.5. High dose hook effect

The assay was tested for a high dose hook effect. Up to a PSA concentration of 2000 ng/mL no hook effect was observed. Please note that if the OD is out of the standard range for highly concentrated samples, the sample must be diluted before the next measurement to obtain correct results.

### 11.6. Correlation

Diametra Total PSA ELISA was compared with the Roche ElecSys total PSA:

$$Y = 0.9644 * X + 0.0741$$

In a second study Diametra Total PSA ELISA was compared to another CE registered PSA ELISA:

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Claims submitted due to customer misinterpretation of laboratory results subject to point 12.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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Ed. 01/2015

DCM138-3

**DiaMetra S.r.l. Headquarter:** Via Calabria 15,  
20090 SEGRATE (MI) Italy

Tel. +39-02-2139184

Fax +39-02-2133354

**Manufactory:** Via Pozzuolo 14

06038 SPELLO (PG) Italy

Tel. +39-0742-24851

Fax +39-0742-316197

E-mail: [info@diametra.com](mailto:info@diametra.com)

## Free PSA

for routine analysis

Enzyme immunoassay for the quantitative determination of free PSA (Prostate Specific Antigen) in human serum and plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 test

REF DKO138

### INTENDED USE

Diametra Free PSA ELISA is used for the quantitative determination of free Prostate Specific Antigen (f-PSA) in human serum or plasma samples. The determination of f-PSA levels is generally used in conjunction with a total PSA (t-PSA) measurement to determine the ratio between f-PSA and t-PSA. This ratio helps to estimate the risk for prostate cancer and to discriminate between elevated t-PSA levels caused by cancerous or non-cancerous conditions. F-PSA determinations are especially recommended for men with elevated t-PSA levels and negative results with digital rectal examination (DRE) in order to decide if a second prostate biopsy is indicated.

### 1. CLINICAL SIGNIFICANCE

Prostate cancer is the most frequent type of cancer found in man and is the second cause of death due to cancer in males. Until recently, digital rectal examination (DRE) was frequently used as only diagnostic modality for the detection of early stages of prostate cancer. In the recent years the determination of serum PSA levels has become the most accepted method to improve the diagnostic specificity of DRE. Although PSA is a tissue specific protein and is not solely tumor specific, it has become the most important marker for prostate carcinoma, showing a better specificity than other biochemical markers used in this context (PAP, total alkaline phosphatase, carcinoembryonic antigen, etc.)

In 1979, Wang et al. isolated a specific antigen for normal prostate tissue and called this protein PSA. PSA is a 33 kDa serine proteinase. Immunohistological studies have shown that PSA is localized in the cytoplasm of prostate acinar cells, ductal epithelium and in the secretion on the ductal lumina, present in normal, benign hyperplastic and malignant prostate tissues as well metastatic prostate cancer and in seminal plasma. If the structural integrity of the prostate is disturbed and/or the gland size is increased, the amount of PSA in the blood plasma may become elevated. In the blood plasma, most of the PSA forms complexes with various proteinase inhibitors. Only a small fraction of PSA circulates as free inactive PSA. Basically three major forms of PSA can be distinguished, only two of which are immunoreactive. The predominant form of PSA is a complex with  $\alpha$ 1-antichymotrypsin (ACT-PSA). Inactive free PSA (f-PSA) represents around 10-40% of the immunologically detectable PSA. The total amount of immunoreactive PSA is known as total PSA (t-PSA). PSA complexed with  $\alpha$ -2-macroglobulin cannot be detected by immunological assays and is therefore frequently called occult PSA (o-PSA).

Current methods of screening men for prostate cancer utilize the detection of t-PSA. Levels of 4.0 ng/ml or higher

are strong indicators of the possibility of prostatic cancer and are an indication for follow-up examinations of the patient. However, elevated serum PSA levels are frequently also attributed to benign prostatic hyperplasia, leading to a high percentage of false positive screening results. A potential solution to this problem involves the determination of free PSA levels. Studies have suggested that the percentage of free PSA is lower in patients with prostate cancer than those with benign prostatic hyperplasia. Thus, the measurement of free serum PSA in conjunction with total PSA, can improve specificity of prostate cancer screening in selected men with elevated total serum PSA levels, which would subsequently reduce unnecessary prostate biopsies with minimal effects on cancer detection rates.

### 2. PRINCIPLE

This f-PSA ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA). The microtiter wells are coated with an anti-f-PSA monoclonal antibody, directed towards an epitope of an antigen molecule. An aliquot of patient serum is incubated in the coated well with enzyme conjugated second antibody (E-Ab), directed towards a different region of the antigen molecule. After incubation the unbound E-Ab is washed off. The amount of bound E-Ab is proportional to the concentration of antigen in the sample. After adding the substrate solution, the intensity of colour developed is proportional to the antigen concentration in the sample. The measured ODs of the Calibrators are used to construct a calibration curve against which the unknown samples are calculated.

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#### 3.1. Reagents and materials supplied in the kit

1. PSA Calibrators (6 vials, 0.5 mL each)  
CAL0 / Sample Diluent (10 mL) **REF DCE002/13806-0**  
CAL1 **REF DCE002/13807-0**  
CAL2 **REF DCE002/13808-0**  
CAL3 **REF DCE002/13809-0**  
CAL4 **REF DCE002/13810-0**  
CAL5 **REF DCE002/13811-0**

Contain non-mercury preservative.

2. Controls (2 vial, 0.5 mL each)  
Low Control **REF DCE045/13804-0**  
High Control **REF DCE045/13805-0**

For concentration see the Certificate of Analysis.  
Contains non-mercury preservative.

3. Coated Microplate (1 breakable microplate)  
Monoclonal antibody anti PSA adsorbed on microplate  
**REF DCE002/13803-0**

4. Assay Reagent (1 vial, 6 mL)  
Contains non-mercury preservative  
**REF DCE043/13843-0**

5. Enzyme Conjugate (1 vial, 6 mL)  
Monoclonal antibody anti PSA conjugated to horseradish peroxidase. Contains non-mercury preservative.  
**REF DCE002/13802-0**

6. 40X Wash Buffer (1 vial, 25 mL)  
Dilute with distilled water before use.  
**REF DCE057/13857-0**

7. TMB Substrate (1 vial, 12 mL)  
Contains TMB (tetramethylbenzidine) and H<sub>2</sub>O<sub>2</sub>  
**REF DCE004/13804-0**

8. Stop Solution (1 vial, 14 mL)  
Contains sulphuric acid (H<sub>2</sub>SO<sub>4</sub>) 0.5M (avoid contact with the Stop Solution, it may cause skin irritations and burns)  
**REF DCE005/13805-0**

#### 3.2. Material required but not provided

- Precision micropipettes (volume: 25 µL and 100 µL) with disposable tips
- Distilled water
- ELISA photometer with 450 nm- and 630 nm-filters
- Timer with 60 min. range or higher
- Microplate washer (optional)
- Vortex or similar mixing tools
- Container for the proper handling of waste and samples after use

#### 4. STORAGE AND STABILITY

- Store the kit and components at +2 °C to + 8 °C
- Bring to room temperature (18-25°C) at least 30 minutes before use. After use put back into the refrigerator. Avoid long time storage at room temperature.
- Do not use the kit or components after the expiry date. For expiry date of the original packed kit see kit label.
- Close the bottles immediately after use.
- Store the plate incl. desiccant in the provided zip-lock pouch. Modules that are not used should always be stored under this condition.
- Ensure that kit components do not freeze.
- Opened kits retain activity for 8 weeks if stored as described above.

### 5. PRECAUTIONS

- ELISA kits are only for in vitro diagnostic use by professionals.
- Serum and plasma samples should be treated as potentially infectious materials. Wear gloves and proper laboratory attire when handling sample materials. Do not eat, drink or smoke in areas where specimen or kit reagents are handled. Do not pipette with the mouth. In case of skin contact, wash with a germicidal soap and copious amounts of water. Seek medical advice if indicated.
- The PSA Calibrators and Controls are of human origin. They have been tested and confirmed negative for HIV, HBsAg and HCV. However, all Calibrators should be treated as potential biohazards.
- Due to the potentially infectious character of samples and kit components all materials that have come in contact with these materials should be sterilized and disposed of according to local regulations. This also includes the liquid waste.
- The assay reagents contain preservatives, TMB, H<sub>2</sub>O<sub>2</sub> or sulphuric acid and may be harmful if ingested. A direct skin or mucosa contact should be avoided. In case of skin contact, wash thoroughly with water and seek medical attention if required.
- The stop solution contains H<sub>2</sub>SO<sub>4</sub>. Since the H<sub>2</sub>SO<sub>4</sub> used to terminate the color reaction is corrosive, the instrumentation employed to dispense it should be thoroughly cleaned after use.
- Do not interchange reagents from different lots or different suppliers.
- Avoid reagent or sample carry-over by using fresh tips for solutions and samples.
- Do not use test kit if zip lock pouch or bottles have been damaged.

### 6. GUIDELINE FOR SAMPLE COLLECTION, PREPARATION AND STORAGE

#### 6.1. Sample collection

Blood samples are collected by veinpuncture. As different factors could influence the PSA level in blood, doctors should ensure that the patient has avoided the following conditions before taking the blood sample.

The following conditions may lead to an increase of PSA levels:

- biking
- sexual intercourse (ejaculation)
- Manipulation of the prostate during medical examinations like DRE, transrectal prostatic ultrasound etc. Please note that especially DRE might lead to an increase of f-PSA only, so that PSA carcinomas might be overlooked
- Prostatitis
- Liver dysfunction

The following conditions may lead to a decrease of PSA levels

- Intake of 5-alpha-reductaseinhibitors, antiandrogens, or GnRH analoga

#### 6.2. Sample preparation

The preparation of serum or plasma samples is performed according to standard techniques. Serum or plasma should be prepared as soon as possible to avoid hemolysis and to improve the stability of PSA.

### 6.3. Storage of samples

For the assay either fresh serum or plasma samples can be used. F-PSA is not as stable as t-PSA.

Literature (Thomas, 2008) recommends to measure samples within 24 hours. In case of longer storage, freeze at -20 °C. A repeated freezing and thawing of samples should be avoided.

#### Note

- Highly lipemic or hemolytic samples can give incorrect analytical results.
- Samples must be free of microbial contaminations.
- Samples containing high titers of rheumatoid factor and human anti-mouse antibodies (HAMA) could give erroneous results.

## 7. ASSAY PROCEDURE

### 7.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are calibrated against WHO 96/668 and have the following concentrations:

	C <sub>0</sub>	C <sub>1</sub>	C <sub>2</sub>	C <sub>3</sub>	C <sub>4</sub>	C <sub>5</sub>
ng/mL	0	0.75	1.50	3	6	12

Contain a preservative.

### 7.2. Preparation of the Wash Solution

Before starting the assay the wash buffer must be diluted to the right concentration. Per well approximately 2 mL of diluted buffer are needed. Calculate the volume of buffer required for the testing. Take 1/40 of the volume wash buffer concentrate and dilute with 39/40 of the volume distilled water.

### 7.1. Test procedure

Note: It is highly recommended to perform all measurements as duplicates. An independent calibration curve should be made for each series of measurements. For best results it is important that the solutions are always added to the wells in the same order to minimize variations in incubation times.

- Prior to use bring all reagents, Calibrators, controls, and samples to room temperature (18-25°C).
- Check that all components are not expired and take care that bottles and plate (inclusive pouch) are not damaged.
- Format the required microplate wells. Keep in mind that all measurements should be performed as duplicate. Document position of wells and respective samples, Calibrators and controls to ensure later identification. Put any unused microwell modules back into the zip lock bag with the desiccant, seal bag and store at 2-8 °C.
- Samples with an expected f-PSA value higher than 12 ng/mL should be diluted with the sample diluent.

Reagent	Calibrator	Sample	Controls
Calibrators C <sub>0</sub> -C <sub>5</sub>	25 µL		
Controls			25 µL
Sample		25 µL	
Assay Reagent	50 µL	50 µL	50 µL
Mix by moving plate on the table (10 sec). Incubate 1 hour at room temperature (18-25°C).			
Enzyme Conjugate	50 µL	50 µL	50 µL
Mix by moving plate on the table (10 sec). Incubate 1 hour at room temperature (18-25°C). Remove solution from the wells by aspirating the liquid or by decanting it. If decanting, tap plate on adsorbent paper to remove residual liquid.  For washing fill plate with diluted wash buffer and wait 10sec before removing the buffer; repeat wash 5 times (for a total of 6 times).  We recommend the following procedure: wash wells 6-times with 250-300 µL/well diluted wash buffer. Preferably use an automated washing procedure. If washing manually, take care that the wash buffer remains in each well for the same time. This is necessary to receive lowest possible CV-values.			
Substrate Solution	100 µL	100 µL	100 µL
Incubate 15 minutes at room temperature (18-25°C) in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read absorbencies (OD) at 450 nm (blinking 630nm)			


### 7.1. Results

1. Calculate the mean absorbance for each duplicate.
2. Subtract the absorbance value of the Calibrator 0 from the mean absorbance values of Calibrators, control and samples.
3. Draw the calibration curve on lin-lin or log-log graph paper by plotting absorbance values of Calibrators against appropriate PSA concentrations.
4. Read off the f-PSA concentrations for the control and the samples.

## 8. VALIDITY OF THE ASSAY

1. The OD 450 nm of the blanking well is lower than 0.150. Higher values indicate a chromogen/substrate contamination. In such a case, repeat the assay carefully checking the reagent.
2. The OD 450 nm of the highest Calibrator (12 ng/mL) must be higher than 0.9. Lower values indicate kit or control decay. In such a case, check the expiry date of the kit before repeating the assay.
3. The control provided should not differ by more than 15% from the concentration stated on the label of the vial if run at least in duplicate.
4. Worksheet and calibration curve of typical assay: not to be used for calculation of actual test results.

Wells	Identity	A 450 nm		Conc. ng/mL
1-2	Cal 0 ng/mL	0.011	0,009	
3-4	Cal 0.75 ng/mL	0.109	0.096	
5-6	Cal 1.5 ng/mL	0.185	0.179	
7-8	Cal 3 ng/mL	0.393	0.392	
9-10	Cal 6 ng/mL	0.871	0.837	
11-12	Cal 12 ng/mL	1.901	1.844	
13-14	Control	0.218	0.231	2.05

Note that the absolute OD values for the Calibrators might vary due to temperature influences or age of the conjugate. As long as the OD values form a calibration curve and remain within the specifications and the control shows the expected value, results for unknown PSA samples are valid.

## 9. QUALITY CONTROL

- It is recommended that internal controls are used in every assay.
- Control results should be within established ranges and should represent preferably low, medium and high concentrations.

## 10. EXPECTED VALUES LIMITATIONS OF THE PROCEDURE

When total PSA values are between 4-10 ng/mL the ratio of f-PSA to t-PSA might be employed to increase the diagnostic specificity of PSA testing. Generally the likeliness of a cancerous condition increases the lower the f-PSA portion is.

A ratio of f-PSA to t-PSA of 25 % usually indicates a high probability of BPH (benign prostate hyperplasia). Low f-PSA likely signals prostate cancer. Most men with prostate cancer have an f- PSA value below 15 %. If free PSA is below 7% prostate cancer is most likely. According to the American Cancer Society and National Cancer Institute men with f-PSA below 7 % should undergo biopsy.

Please note that an isolated f-PSA concentration is of no diagnostic value.

For a detailed discussion of probabilities of prostate cancer and percent free PSA also see Gion et al. (1998).

The above values are only for user's guidance. If possible, it is recommended for each laboratory to establish its own specific values that take into consideration a population indigenous to the area where the laboratory is located.

**Note: PSA values and the f-PSA to t-PSA can only be used to estimate the cancer risk. They should always be interpreted in conjunction with other clinical findings and should not be used as a sole basis for prostate cancer diagnosis.**

## 11. PERFORMANCE CHARACTERISTICS

### 11.1. Detection limit

The limit of detection for this kit is 0.1 ng/mL.

### 11.2. Precision

Intra- and inter-assay precisions were established by analysing three patient sera of different PSA concentrations.

The results are summarized in Tables 1 and 2.

Table 1) Intra-assay precision

Patients	Number of replicates	Mean ng/mL	SD ng/mL	CV %
1	10	0.31	0.024	7.7
2	10	1.75	0.107	6.1
3	10	9.63	0.623	6.5

Table 2) Inter-assay precision

Patients	Number of replicates	Mean ng/mL	SD ng/mL	CV %
1	10	2.86	0.17	6.1
2	10	0.75	0.06	7.8
3	10	2.1	0.16	7.8

### 11.3. Recovery

Spiked serum samples were prepared by adding aliquots of samples with highly elevated free PSA to normal male serum samples. The recovery of the antigen was in between 91.8 – 116% and in average 102%.

Table 3) Recovery

Sample	Expected value (ng/mL)	Observed value (ng/mL)	Recovery %
1	5.48	5.47	99,8
	2.74	3.19	116
	5.48	5.03	91.8
2	2.78	2..81	101
	5.57	5.68	102
3	2.74	2.82	103

### 11.4. Hook effect

No hook effect has been noticed with samples up to 5000 ng/mL.

### 11.5. Correlation

The f-PSA ELISA (EIA-4189) was compared to a commercially available and CE marked two step Free PSA Elisa:

$$Y = 0.8922x + 0.11; R^2 = 0.8977$$

In a second study the f-PSA ELISA (EIA-4189) was compared to a second reference f-PSA ELISA:

$$Y = 1.0708 x; R^2 = 0.7969$$

### 11.6. Specificity

The assay is highly specific for free PSA, with a relatively low cross-reactivity to other proteins and polypeptides, lipids or chemotherapeutic agents that might be present in patient samples.

Table 4) Specificity

Antigens	Amount added	Cross reaction
<b>Proteins</b>		
AFP	10 µg/mL	No
CEA	10 µg/mL	No
HCG	10 µg/mL	No
Lactalbumin	10 µg/mL	No
PAP	1 µg/mL	No
<b>Interfering substances</b>		
Bilirubin	0.2 mg/mL	No
Triglyceride	15 mg/mL	No
Hemoglobin*	0.1 mg/ml	No
<b>Chemotherapeutic Agents</b>		No
Cyclophosphamid	800 µg/mL	No
Doxorubicin * HCl	20 µg/mL	No
Diethylstilbestrol	2 µg/mL	No
Flutamide	10 µg/mL	No
Methotrexate	50 µg/mL	No

\* at higher concentration hemoglobin results in too high OD values, hemolytic samples should thus be avoided.

## 12. LEGAL ASPECTS

### 12.1. Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact Diametra.

### 12.2. Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 12.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 12.3. Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 14.2. are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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