Please read the Instructions for Use carefully before using the product.

**Intended use:**

This product is used for hemodialysis treatment of acute or chronic renal failure. It forms a circuit with bloodlines to complete the treatment. The product directly contacts with patient's blood during treatment. Blood flows into the hollow fiber membrane, and dialysate flows outside of the membrane, forming osmotic pressure on both sides, resulting in the removal of small molecular toxins from the patient, meanwhile forming a transmembrane pressure to remove extra fluid from the patient.

**Intended patient population:**

This product is used in patient undergoing hemodialysis treatments. When selecting this product for a treatment, the total extracorporeal blood volume (i.e. in combination with the dialyzer, the bloodline set and any other accessories) shall not exceed 10% of the patient's blood volume.

**Intended user / environment:**

This product shall be handled by healthcare professionals.

Application of aseptic techniques shall be ensured throughout the entire procedure.

**Contraindications:** Special contraindications are unknown. Contraindications for dialysis are applicable. Strict monitoring shall be performed for the patients who have tendencies of blood bleeding or clotting during the treatment. Treatment shall be terminated in the event of any complications affecting patient's condition.

**Side effects:** Adverse reactions such as hypertension, hypotension, headache, and nausea may be associated with hypovolemia or hypervolemia and can usually be avoided with management of patient fluid, electrolyte balance, blood flow rate, and ultrafiltration rate.

In rare cases, thrombocytopenia or hypersensitivity reactions including anaphylactic or anaphylactoid reactions to the dialyzer, or other elements in the extracorporeal circuit, may occur during hemodialysis. Hypersensitivity reactions may cause mild to severe signs and symptoms, including: itching, flushing, hives, swelling, fever, leukopenia, hypotension, hypertension, shortness of breath with wheezing, arrhythmias, and/or respiratory arrest.

Patients with a history of hypersensitivity reactions or patients who have a history of being highly sensitive and allergic to a variety of substances should be carefully monitored during treatment. With severe hypersensitivity reactions, dialysis must be discontinued and aggressive first line therapy for hypersensitivity reactions must be initiated. The decision to return the patient's blood in the event of a hypersensitivity reaction is determined by the physician.

**Cautions and Warnings:**

- The dialysate and the blood should flow in counter-current.
- The maximum flow rate of dialysate is 800mL/min.
- Not recommended for paediatric use and **DO NOT** use on non-de-aerated dialysis fluid delivery systems.
- During the treatment, the transmembrane pressure shall not exceed 66.66kPa (500mmHg).
- The blood flow rate shall not be lower than 150mL/min, but not higher than 500mL/min.
- Use an aseptic technique when preparing the circuit and the dialyzer for use to avoid contamination.
- Hollow Fiber Dialyzer shall be used under medical supervision or adequately trained personnel. It is suggested to use together with hemodialysis machine, dialysis water, concentrated solution and dialysate which are in accordance with the international standards to avoid bacterial and pyrogenic contamination.
- The fluid pathways (blood and dialysate) is sterile and non-pyrogenic, sterilized by irradiation. **DO NOT** use the expired product. **DO NOT** use the device if package is damaged or if protective caps are not in place.
- This dialyzer is for single use and reuse is strictly prohibited. Reprocessing of this product may lead to adverse patient reactions and/or device failure. It should be discarded according to local requirements relevant to disposal of infectious medical waste to prevent contamination.
- If abnormal conditions arise during the dialysis, such as bubbles, foreign matter, blood leak, or clotting, proper measures shall be taken according to doctor's advice.
- Please use only with ultrafiltration controlling dialysis system.
- No use of dialysis supply systems without degassing.
- Do not use expired products. Please check the expiration date on label before use.
- Ensure the connection between bloodline and dialyzer is strictly tight. Do not use this product if the bloodline connectors cannot fit for the dialyzer. Make sure all the connectors are tight during treatment to prevent blood leakage or any air entry
- If the device is used below a blood flow rate of 150mL/min, or below a dialysate flow rate of 500mL/min, it may have diminished performance.
- If serious incident occurs, please inform the manufacturer or local competent authority.
- Please stop using if any component causes allergic reaction.

**Instructions for Installation**

- Take out the dialyzer from the pouch, and check whether the dialyzer and its components are in good condition.
- Set up the dialyzer vertically to the holder.
- Make sure the dialyzer stays firmly in the holder.

**Recommended usage method**

**1. Priming**

- Prepare no less than 500mL normal saline or dialysate and add appropriate amount of heparin if needed under physician's prescription.
- Place the arterial line and the venous line onto the dialysis machine according to the instructions for use of the extracorporeal blood circuit.
- Connect the arterial line, the venous line and the dialyzer.
- Control the flow rate of blood pump within 80-100mL/min, use normal saline or dialysis fluid until totally remove air from blood line and blood compartment of the dialyzer. The direction of normal saline or online dialysis fluid flow is arterial line→dialyzer→venous line, and counter-current is forbidden.
- Turn the flow rate up to 200-300mL/min, connect the dialysate connector to the dialyzer. Turn the dialyzer so that the venous end is downward, and make sure to remove all air from the dialysate compartment.

**Note:** Check for the presence of air in the blood circuit before patient connection and run extra priming/ recirculation/ ultrafiltration if any air may remain.

**2. Anticoagulation**


- Carry out the anticoagulation according to prescription.

**3. Termination of Treatment**

- For reinfusion and completion of the treatment, follow the instruction on the dialysis machine in use.
- Turn the flow rate up to 50-100mL/min, use adequate saline or substitution fluid to complete blood reinfusion. Turn off the blood pump and disconnect from patient's venous line.


**Note:** **DO NOT** turn off the air monitoring system before return blood to the patient completely to prevent the air from entering into the patient from the blood line.

**A typical connection diagram**

**Main Structure:**

	A	B	C	D	E
<b>Component</b>	Housing	Dialysis Membrane	Cover	Sealant	Protective Cap
<b>Material</b>	Polypropylene	PES Membrane	Polypropylene	Polyurethane	Polyethylene

	Do not re-use		A sterile fluid path that has been sterilized using irradiation		Fluid path is non-pyrogenic		CE marking
	Date of manufacture		Manufacturer		Temperature limit		Medical device
	Batch code		Use-by date		This end up		Single sterile barrier system
	Catalogue number		Keep dry		Fragile, handle with care		Single sterile barrier system with protective packaging outside
	Do not use if package is damaged and consult instructions for use		Keep away from sunlight		Consult instructions for use		Unique Device Identifier
	Do not re-sterilize		Humidity limitation		Authorized representative in the European Community/European Union		

**Product Performance:** This dialyzer has reliable performance, which can be used for hemodialysis. The basic parameters of product performance and the laboratory data of this series will be provided as follows for reference.

**Note:** The laboratory data of this dialyzer was measured according to the standards ISO 8637-1.

These data represent typical *in vitro* performance. *In vivo* performance will differ due to the patient's blood composition and clinical settings.

Model	B-08PF			B-09PF			B-10PF			B-11PF			B-12PF			B-13PF			B-14PF			B-15PF			B-16PF					
Product ref.	B-08PF			B-09PF			B-10PF			B-11PF			B-12PF			B-13PF			B-14PF			B-15PF			B-16PF					
Test Conditions: Q <sub>0</sub> = 500 mL/min, Temperature: 37°C±1°C, Q <sub>1</sub> = 10 mL/min																														
Clearance / Q <sub>0</sub> (mL/min)	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400
Urea (mL/min)	182	227	256	183	231	260	184	235	264	185	239	269	187	243	273	188	247	277	189	251	281	190	255	285	191	259	290			
Creatinine (mL/min)	168	197	236	170	203	240	172	209	244	174	215	248	176	221	252	178	226	256	180	232	260	181	238	267	183	244	275			
Phosphate (mL/min)	161	193	232	164	199	238	167	205	244	170	211	250	173	216	256	176	222	261	179	228	267	182	234	273	185	240	280			
Vitamin B <sub>12</sub> (mL/min)	90	88	101	95	96	106	99	104	113	103	111	121	108	118	128	111	125	135	114	132	142	118	140	150	123	148	158			
Pressure drop of blood compartment (mmHg)	<55	<75	<95	<55	<75	<95	<55	<75	<95	<55	<75	<95	<55	<75	<95	<50	<70	<90	<50	<70	<90	<50	<70	<90	<45	<65	<85			
KUF (mL/hr/mmHg) Q <sub>0</sub> =300 mL/min	19			19			20			21			22			23			24			25			26					
Priming volume (mL)	49			55			61			67			73			79			85			91			98					
Effective membrane area (m <sup>2</sup> ).	0.8			0.9			1.0			1.1			1.2			1.3			1.4			1.5			1.6					

Model	B-17PF			B-18PF			B-19PF			B-20PF			B-21PF			B-22PF			B-23PF			B-24PF			/					
Product ref.	B-17PF			B-18PF			B-19PF			B-20PF			B-21PF			B-22PF			B-23PF			B-24PF			/					
Test Conditions: Q <sub>0</sub> = 500 mL/min, Temperature: 37°C±1°C, Q <sub>1</sub> = 10 mL/min																														
Clearance / Q <sub>0</sub> (mL/min)	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	200	300	400	/	/	/
Urea (mL/min)	192	262	295	193	266	300	194	270	305	195	274	311	196	279	316	197	284	321	197	287	325	198	291	330	198	291	330	/	/	/
Creatinine (mL/min)	184	250	281	186	256	287	187	261	293	189	267	299	190	273	305	192	280	311	193	283	316	195	286	322	195	286	322	/	/	/
Phosphate (mL/min)	187	245	286	190	250	292	191	255	298	192	260	305	193	266	311	195	272	317	196	276	323	197	281	329	197	281	329	/	/	/
Vitamin B <sub>12</sub> (mL/min)	126	155	166	130	162	175	134	169	181	138	177	188	142	183	195	146	190	202	149	196	208	153	201	214	153	201	214	/	/	/
Pressure drop of blood compartment (mmHg)	<45	<65	<85	<40	<60	<80	<40	<60	<80	<40	<55	<75	<40	<55	<75	<40	<55	<75	<40	<55	<75	<40	<55	<75	<40	<55	<75	/	/	/
KUF (mL/hr/mmHg) Q <sub>0</sub> =300 mL/min	27			28			29			30			31			32			33			34								
Priming volume (mL)	104			110			117			125			131			137			144			150								
Effective membrane area (m <sup>2</sup> ).	1.7			1.8			1.9			2.0			2.1			2.2			2.3			2.4								

**Transportation and storage:** Handle the product with care during transportation. Avoid direct exposure to sunlight, rain or snow. Store this product in a well-ventilated indoor space without corrosive gases. Store at temperature of 0°C–40°C and not more than 80% relative humidity. DO NOT store this product together with chemicals and/or moist articles.

**Shelf life:** please refer to the actual labeling.

**After-sale service:** Please keep the original packaging for any investigation on product quality.

**EU REP** <EU Representative>

**MT Promed Consulting GmbH**  
 Add: Ernst-Heckel-Straße 7, 66386 St. Ingbert, Germany  
 Tel: +49 (0) 6894 581020 Fax: +49 (0) 6894 581021

<Manufacturer>

**Bain Medical Equipment (Guangzhou) Co., Ltd.**  
 Add.: No.10, Juncheng Road, Eastern Area, Economic and Technological Development District, 510760 Guangzhou, P.R. China.  
 Tel: +86-20-66856868 Fax: +86-20-32067500  
 E-mail: [sales@baingz.com](mailto:sales@baingz.com)

Keep this Instructions for Use until all of the products in this carton have been used up.

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