Generic Name : Magnesium Kit (Xylidyl Blue Method)**Abbreviated name : Mg****Order Information**

Cat. No.	Package size
MG0102	R 4×40 mL
MG0103	R 6×40 mL
MG0104	R 4×45 mL
MG1102	R 1×25 mL
MG1103	R 1×25 mL
MG1104	R 1×25 mL
MG0105	R 4×250 mL

Intended use

In vitro test for the quantitative determination of Mg concentration in serum, plasma and urine on photometric systems.

Summary ¹

Magnesium is one of the most abundant cations in the body, and plays an important role in cell respiration, glucose metabolism and transmembrane transport.

Magnesium can activate more than 300 kinds of enzymes, the famous of which is Na⁺-K⁺-ATP enzyme.

Hypomagnesemia may cause by defective gastrointestinal absorption, body fluid losses, diuretic therapy, aminoglycoside treatment, endocrinopathy and alcoholism, and the inherited disease is also an infrequent reason.

Hypermagnesemia is found in acute and chronic renal failure, magnesium excess, dehydration and diabetic acidosis.

Method

Xylidyl blue method

Reaction Principle

Xylidyl blue + Magnesium $\xrightleftharpoons{\text{OH}^-}$ Xylidyl blue-Magnesium complex

By using the EGTA to eliminate the interference of calcium, magnesium ions combine with xylidyl blue to produce a xylidyl blue -magnesium complex at an alkaline solution. The absorbency increase is directly proportional to the concentration of magnesium.

Reagents**Components and concentrations**

R_i	Ethanolamine	49 mmol/L
	EGTA	0.13 mmol/L
	Xylidyl blue	0.09 mmol/L
	Surfactant	<2% (m/v)

Warnings and precautions

1. For in vitro diagnostic use only.
2. Take the necessary precautions for the use of laboratory reagents.
3. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
4. Disposal of all waste material should be in accordance with local guidelines.
5. Material safety data sheet is available for professional user on request.

Reagent Preparation

Single reagent is ready to use.

Storage and stability

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

Once opened, the reagents are stable for 10 days when refrigerated on the analyzer or refrigerator.

Contamination of the reagents must be avoided.

Do not freeze the reagents.

Reagent blank absorbency

The absorbance of reagent blank at 546 nm should be <1.4 A.

Materials required but not provided

1. Calibrator and controls as indicated below.
2. NaCl solution 9 g/L.
3. General laboratory equipments.

Specimen collection and preparation ²

1. Serum, heparin plasma and urine are suitable for samples. Whole blood, hemolysis not recommended for use as a sample. Freshly drawn serum is the preferred specimen.
2. Collect urine sample in a metal-free container with no preservatives; bring the urine to pH 3-4 prior to assay.
3. Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.
4. Centrifuge samples containing precipitate before performing the assay.
5. Stability:

Plasma must be assayed fresh.

Serum: 7 days at 20-25°C

7 days at 4-8°C

1 year at -20°C

Urine: 7 days at 20-25°C

7 days at 4-8°C

1 year at -20°C

Assay procedure

	Blank	Sample
R	1000 µL	1000 µL
Dist. water	10 µL	—
Sample	—	10 µL

Mix thoroughly at 37°C, and read the absorbance 5 min. later.

$$\Delta A = [\Delta A \text{ sample}] - [\Delta A \text{ blank}]$$

Application sheets for BS series analyzers are available in this document. Refer to the appropriate operator manual for the analyzer-specific assay instructions.

Calibration

1. It is recommended to use the Human multi-calibrator from Mindray and 9 g/L NaCl for two-point calibration. Traceability of the multi-calibrator can refer to the calibrator instructions for use of Mindray Company.

2. Calibration frequency:

If reagent bottle is onboard the analyzer for more than 5 days

After reagent lot changed.

As required following quality control procedures.

Quality control

At least two levels of control material should be analyzed with each batch of samples. In addition, these controls should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual. We recommend using the Human Assayed Control made by Mindray to verify the performance of the measurement procedure; other suitable control material can be used in addition.

Each laboratory should establish its own internal quality control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

Calculation

The analyzer calculates the concentration of each sample automatically after calibration.

Conversion factor: $\text{mg/dL} \times 0.4113 = \text{mmol/L}$

Or: $C \text{ sample} = (\Delta A \text{ sample} / \Delta A \text{ calibration}) \times C \text{ calibration}$

Reference Intervals ¹

Each laboratory should establish its own reference intervals based upon its patient population. The reference intervals measured at 37°C listed below were taken from literature:

Sample Type		S.I. Units
Serum / Plasma	Neonates	0.48–1.05 mmol/l
	Children	0.60–0.95 mmol/l

	Women	0.77–1.03 mmol/l
	Men	0.73–1.06 mmol/l
Urine	Male	3–5 mmol/24 h

Performance Characteristics

Representative performance data obtained from Mindray system (Mindray BS series analyzers / Mindray Mg Reagent) is given below. Results may vary if a different instrument, an individual laboratory or a manual procedure is used.

Limitations-interference

The following substances were tested for interference with this methodology. Criterion: Recovery within ± 10 % of initial value.

Substance	Level Tested	Observed Effect
Ascorbic acid	30 mg/dL	NSI*
Bilirubin	40 mg/dL	NSI
Lipemia	200 mg/dL	NSI
Hemoglobin	250 mg/dL	NSI

* NSI: No Significant Interference (within ± 10 %)

Linearity range

The Mindray System provides the following linearity range:

Sample Type	Conventional Units	S.I. Units
Serum / Plasma / Urine	0.10-4.16 mg/dL	0.04–2.05 mmol/L

If the value of sample exceeds 2.05 mmol/L, the sample should be diluted with 9 g/L NaCl solution (e.g. 1+ 1) and rerun; the result should be multiplied by 2.

Analytic Sensitivity/Limit of Detection

The lowest measurable Mg concentration that can be distinguished from zero is 0.04 mmol/L with 99.7% confidence.

Precision

Precision performance using the CLSI Approved Guideline EP5-A2 to assay serum control appears in the table below³.U: mmol/L

Type of Imprecision	Level I			Level II		
	Mean	SD	CV %	Mean	SD	CV %
Within-run	0.867	0.005	0.550	1.370	0.006	0.427
Between-run		0.012	1.347		0.017	1.256
Between-day		0.010	1.171		0.015	1.078
Within-device		0.016	1.867		0.023	1.709

Method Comparison

A comparison between Mindray System (Mindray BS series analyzers /Mindray Mg Reagent) (y) and Hitachi/Roche System (Hitachi /Roche Mg) (x) using 61 samples gave following correlation (mmol/L): $y=1.012x+0.0225$,

$R^2 = 0.9982$

Details of the comparison experiments are available on request.

References

1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998.
2. Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag, 1996.
3. CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

Graphical symbols

In Vitro Diagnostic
medical device

Batch Code

European
Conformity

Authorized representative in
the European Community

Use By

Consult
Instructions for use

Temperature
Limit

Manufacturer

Catalogue
number

© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R.China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp.
GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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Referências

1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998.
2. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the Patient to the Laboratory. Darmstadt: GIT Verlag, 1996.
3. CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P.R.China

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© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tüm hakları saklıdır
Üretici:Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
AT Temsilcisi:Shanghai International Holding Corp.GmbH(Avrupa)