Medtronic

Amplatz Goose Neck™

Snare Kit
Комплект на примка
Souprava zachytávacího očka
Snareset
Schlingenset
Κιτ βρόχου
Kit de lazo
Lingukomplekt
Slingasarja
Kit de laso
Komplet omče
Elfogőeszköz-készlet
Kit Cerut
Kit del dispositivo di recupero a cappio
Ілмек жинағы
홀가미 키트
Gaudyklės rinkinys
Satvērēja komplekts
Комплет со јамка
Teugelset
Snaresett
Zestaw pętli
Conjunto de ansa
Trusă cu ansă
Набор ловушки
Súprava zachytávacieho očka
Komplet zanke
Komplet sa zamkom
Snarsats
Kement Seti
Набір пастки
Bộ dụng cụ thòng lọng

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Petunjuk Penggunaan • Istruzioni per l'uso • Қолдану жөніндегі нұсқаулар • 사용 지침 • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з використання • Hướng dẫn sử dụng

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English

1 Device description

The Amplatz Goose Neck snare kit contains 1 snare, 1 snare catheter, 1 snare introducer, and 1 torque device. The snare is made of a nitinol cable and a gold-plated tungsten loop. The preformed snare loop can be introduced through catheters without risk of snare deformation because of the snare's superelastic construction. The snare catheter contains a platinum-iridium, radiopaque marker band.

2 Intended purpose

The intended purpose of the Amplatz Goose Neck snare kit is to retrieve and manipulate foreign bodies in the cardiovascular system during interventional and manipulation procedures.

2.1 Intended patient population

The Amplatz Goose Neck snare kit is intended for use in adult patients requiring retrieval and manipulation of foreign bodies located in the cardiovascular system.

2.2 Indications for use

The Amplatz Goose Neck snare kit is intended for use in the cardiovascular system to retrieve and manipulate foreign bodies.

2.3 Contraindications

The Amplatz Goose Neck snare is contraindicated for the following scenarios:

- Removing foreign bodies that have become entrapped by tissue growth.
- Fibrin sheath stripping in the presence of atrial or ventricular septal defects or persistent foramen ovale.
- Removal of implanted pacing leads.

2.4 Clinical benefits

The device provides clinical benefit by assisting in the retrieval or manipulation of lost, broken, or migrated intravascular objects, which minimizes various vascular and systemic complications.

2.5 Intended users

The Amplatz Goose Neck snare kit should only be used by trained physicians who have experience with interventional and manipulation procedures in the cardiovascular system.

The device should only be used within a clinical setting in a sterile environment.

2.6 Performance characteristics

The snare is a nitinol cable with a single, gold-plated tungsten loop. The 90° angle design of the loop allows for wall apposition and the retrieval and manipulation of foreign bodies. The gold-plated tungsten loop allows for visibility under fluoroscopy. The nitinol in the device provides a superelastic construction, allowing introduction of the snare through a catheter without the risk of snare deformation.

3 Warnings

- This device is not intended for the removal of foreign bodies entrapped by tissue growth.
- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.

- Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.
- The device is made of materials that contain nickel. Some patients are allergic to, or become sensitive to, nickel.

5 Potential adverse events

The potential adverse events (or complications) that may occur or require intervention with the use of this device include, but are not limited to, the following conditions:

- Allergic reaction
- Arrhythmia
- Cardiac injury
- Complications due to damage to an indwelling catheter or other foreign body
- Device entrapment
- Embolism
- Infection
- Myocardial infarction
- Prolonged procedure
- Pulmonary embolism
- Stroke
- Vessel damage
- Vessel perforation

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

After the EUDAMED website is launched, the Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 0763000B00003667T.

6 How supplied

The device is supplied sterile and is sterilized using ethylene oxide.

7 Storage

Store the device in a dry area away from sunlight.

8 Directions for use

Warning: This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.

8.1 Preparation

1. Select the appropriate size snare for the location of the foreign body. In general, the snare loop should approximate the size of the vessel in which the snare will be used.
2. Remove the snare and snare catheter from their protective holders and inspect for any damage.
3. Flush the guidewire lumen with sterile heparinized saline.

8.2 To insert the snare and snare catheter together as a unit

1. Remove the introducer and torque device from the proximal end of the snare shaft.
2. Load the snare into the snare catheter by inserting the proximal (nonlooped) end of the snare into the distal (nonhubbed) end of the snare catheter until the proximal end of the snare shaft exits the hub and the loop end can be retracted into the distal end of the snare catheter.
3. Inspect the device by extending and retracting the snare loop through the distal end of the catheter 2 or 3 times, while carefully examining the catheter, the radiopaque marker band, and the snare for defects or damage.
4. Use a guiding catheter or guidewire to insert the snare and snare catheter together.

8.2.1 With an indwelling guide catheter

Pull the snare just within the distal end of the snare catheter and advance both together through the indwelling guide catheter until the distal end of the snare catheter is positioned as desired.

8.2.2 With a guidewire

1. Place the loop of the snare over the proximal end of the wire and cinch the snare loop over the wire by advancing the snare catheter over the snare shaft.
2. Advance the snare and the snare catheter together over the wire and into the guide catheter until the distal end of the snare catheter is positioned as desired.

8.3 To insert the snare after the catheter is within the vasculature

1. Move the snare introducer (provided on the proximal end of the snare) distally until the loop of the snare is enclosed within the tubing portion of the introducer.
2. Insert the distal end of the introducer into the hub of the indwelling snare catheter until resistance is felt. Resistance indicates that the tip of the introducer is properly aligned with the inner lumen of the catheter.
3. While keeping the introducer as straight as possible, grasp the shaft of the snare just behind the hub of the introducer and advance the snare until it is secure within the lumen of the catheter. To remove the introducer, first remove the torque device from the snare shaft and then pull the introducer off the proximal end of the snare shaft.

8.4 Procedure

8.4.1 Snare assisted manipulation or retrieval

1. Remove the indwelling balloon or delivery catheter, if present. It may be necessary to change or extend the indwelling guidewire to facilitate balloon removal and to upsize the indwelling guide catheter to accommodate the snare.
2. Position the snare just proximal to the foreign body with one of the methods described in *Section 8.2* or *Section 8.3*.
3. Gently push the snare shaft forward within the snare catheter to move the snare forward and completely open the loop. Advance the loop forward, around the proximal end of the foreign body.
4. While maintaining the position of the snare shaft, advance the snare catheter to cinch the snare loop around the foreign body.
Note: Attempting to close the loop by pulling the snare shaft within the snare catheter moves the loop from its position around the foreign body.
5. To manipulate a foreign body, maintain tension of the snare catheter on the snare loop to retain the hold on the foreign body, and move the snare and snare catheter together to manipulate the foreign body to the desired position.
6. To retrieve a foreign body, maintain tension of the snare catheter on the snare loop and move the snare and snare catheter together proximally to or into the guide catheter. Then withdraw the foreign body through or together with the guiding catheter. To withdraw large foreign bodies, if necessary insert larger guiding catheters or perform a cut down at the peripheral site.

9 Disposal

Caution: Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

10 Disclaimer of warranty

Important: This disclaimer of warranty does not apply in any countries where such a disclaimer is not permitted by law.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic, Inc.
Medtronic B.V.
Medtronic Limited