


**I. Intended Use / Indication**

The Resistant self-expanding peripheral stent system is intended to improve luminal diameter in the treatment of symptomatic *de novo* or restenotic lesions up to 190 mm in length in the native superficial femoral artery and/or proximal popliteal arteries with reference vessel diameters ranging from 4.0 mm – 7.0 mm.

The Resistant self-expanding peripheral stent system is indicated for improving luminal diameter in patients with atherosclerotic disease of the common and/or external iliac arteries up to and including 70 mm in length, with a reference vessel diameter of 8.0 mm – 11.0 mm.

**II. Device Description**

The "RESISTANT" is a self-expanding peripheral Nitinol stent system for permanent implantation. The self-expanding stent is made of a nickel titanium alloy (Nitinol) and is pre-mounted on a 6F, 0,035" over-the-wire delivery system. The delivery system is sterilized with ethylene oxide and comprises a detachable outer tube (5); an inner shaft (2); a self-expanding stent with radiopaque markers on both ends of the stent, which is constrained during application of the delivery system between the inner shaft and the outer tube; a retainer tube (3), to support the stent during deployment; a Y-connector (6) and a wing (10) for deployment of the stent, whereby via the luer-connection (7) at the Y-connector the flushing of the space in-between and via the luer-connection (11) at the wing the flushing of the guide wire lumen is done. A system lock (8) prevents the possibility of a premature deployment. At the distal end of the retainer tube and outer tube radiopaque markers (3 & 4) are located for positioning of the stent.




Selection of stent size is based on angiograms to determine the Reference Vessel Diameter and the Lesion Length. It is necessary to select a stent that has an unconstrained diameter at least 0.5 mm larger than the largest reference vessel diameter (RVD) to achieve secure placement according to the following stent size selection table. A higher difference between RVD and unconstrained stent diameter is leading to a higher cronical outward force towards the vessel.


Unconstrained Stent Ø [mm]	Nominal Stent Length [mm]	Sheath Compatibility	Length Foreshortening	Reference Vessel Diameter (RVD) [mm]
5	20, 30, 40, 60, 80, 100, 120, 150, 200	6F	1.6 %	4.0 – 4.5
6	20, 30, 40, 60, 80, 100, 120, 150, 200	6F	1.8 %	5.0 – 5.5
7	20, 30, 40, 60, 80, 100, 120, 150, 200	6F	2.0 %	6.0 – 6.5
8	20, 30, 40, 60, 80, 100, 120, 150, 200	6F	2.3 %	7.0 – 7.5
9	40, 60, 80	6F	2.4 %	8.0 – 8.5
10	40, 60, 80	6F	2.6 %	9.0 – 9.5
12	40, 60, 80	6F	3.1 %	11.0 – 11.5

**III. Technical Specifications**

**IIIa. Catheter Specifications**




<b>Catheter Design</b>	<b>Over-The-Wire (OTW), Pull Back System</b>
<b>Recommended Introducer Sheath</b>	<b>6F</b>
<b>Maximum Guidewire</b>	<b>0.035"</b>
<b>Useable Shaft Lengths</b>	<b>85 cm / 135 cm (850 mm / 1350 mm)</b>
<b>Tip</b>	<b>Soft conical tip</b>
<b>X-Ray Marker</b>	<b>Distal and Proximal to Stent</b>
<b>Sterilization Method</b>	<b>Ethylen Oxide (ETO)</b>
<b>Shelf Life</b>	<b>3 years</b>



- |                      |                     |
|----------------------|---------------------|
| 1 Soft Tip           | 7 Luer (Flushing)   |
| 2 Inner Tube         | 8 System Lock       |
| 3 X-Ray Marker       | 9 Stiffening Tube   |
| 4 X-Ray Marker       | 10 Hub              |
| 5 Outer Braided Tube | 11 Luer (Guidewire) |
| 6 Y-Connector        |                     |

**IIIb. Stent Specifications**






Ø4 – 7 mm		Ø8 mm	Ø9 – 12 mm	
15 Zigzag elements 3 interlinks		16 Zigzag elements 4 interlinks	12 Zigzag elements 4 interlinks	
High flexibility moderate radial strength		High flexibility moderate radial strength	High radial strength moderate flexibility	
Clinical indication				
Femoral Popliteal Tibial	Bile duct Iliac	Femoral	Bile duct Iliac	Iliac

Design	Open Cell Design		
Material	Nitinol		
Stent Diameter	5 – 12 mm		
Dimensions depending on diameter	5 - 7 mm	8 mm	9 - 12 mm
Strut Thickness	215 µm	215 µm	215 µm
Strut Width (Main Segment)	113 µm	135 µm	160 µm
Strut Width (Interlink)	111 µm	113 µm	220 µm
Surface finish	Electropolishing		
Vessel Coverage	10 % - 22 %		
X-Ray Marker	3 / 4 Tantal marker distal and proximal		

#### IV. Biological Characteristics







The RESISTANT stent system is compliant to EN ISO 10993-1:2010

#### V. Sterilization Characteristics, Shelf Life

	The device is sterilized by EO gas.
	The device has a <b>shelf life of 3 years</b> , if stored under condition specified in clause VI.
	Do not re-sterilize the device

#### VI. Packaging, Transportation and Storage Characteristics

Single sterile, packed in Blister sealed with Tyveck® as primary packaging. The blister is placed together with IFU in a carton box as secondary packaging. The carton box is sealed with labels.

	Do not use the device if the package is damaged
	Do not re-use the device
	Lower and upper temperature limitation
	Keep away from sunlight
	Keep dry
	Fragile, handle with care

#### VII. Material Specification

##### a) Silicone

The device does not contain any silicone.

##### b) Phthalates

The device does not contain any phthalates.

##### c) Latex

The device does not contain any latex.

##### d) PVC

The device does not contain any PVC.


##### e) Animal Tissue

The device does not contain material of animal origin or derivatives of animal origin.

##### f) Human Blood Derivates

The device does not contain any derivatives based on human blood.


VIII. Art Work Packaging Design and Label




Front of the product packaging



Back of the product packaging with label



Back of the blister packaging with label



Front of sterile blister packaging

**IX. Label (Size 140 mm x 210 mm)**

**OTW Peripheral Self-expandable Nitinol Stent System**

**RESISTANT 8 x 80 - 135** LOT 5119352 2018-09

diameter length shaft length

---

<p>Stent diameter Stentdurchmesser Diamètre du stent Diámetro del stent Diámetro do stent Stent çapı диаметр стента</p> <p><b>8 mm</b></p>	<p>Shaft diameter Schaftdurchmesser Diamètre cathéter Diámetro del cuerpo Diámetro da bainha Çevre çapı диаметр шприца</p> <p><b>6 F</b> (2.00 mm)</p>
<p>Stent length Stentlänge Longueur du stent Lunghezza dello stent Lungitud del stent Comprimento do stent Stent uzunluğu длина стента</p> <p><b>80 mm</b></p>	<p>Max. rec. guide wire Max. Führungsdrahtdurchmesser Guide max. rec. Diámetro máximo filo guía Diámetro máx. del hilo conductor Diámetro máximo do fio guia Макс. рекомендованный проводник</p> <p><b>0.035"</b> (0.89 mm)</p>
<p>Catheter length Katheterlänge Longueur du cathéter Lunghezza del catetere Lungitud del catetere Comprimento do catetere Katheter uzunluğu длина катетера</p> <p><b>135 cm</b></p>	<p>Content Inhalt Contenu Contenido Conteúdo İçerik количество</p> <p><b>1</b></p>

+EUCARSCO0801359/33015182455119352B

+EUCARSCO08013556

**LOT** 5119352

**2015-10**

**2018-09**

**REF** RSO080135

**LOT**

**eucatech AG - Rebgartenweg 27  
DE-79576 Weil am Rhein - Germany**

**STERILE EO**

**eucatech AG  
RESISTANT**

**LOT** 5119352

**REF** RSO080135

**eucatech AG  
RESISTANT**

**LOT** 5119352

**REF** RSO080135

**LOT** 5119352

**LOT** 5119352

**REF** RSO080135

**LOT** 5119352

**REF** RSO080135

## X. Size Matrix / Order Information

### ORDER INFORMATION

Shaft length: 85 cm

Stent I.D.	Stent length							
	20mm	40mm	60mm	80mm	100mm	120mm	150mm	200mm
5 mm	RSO05020085	RSO05040085	RSO05060085	RSO05080085	RSO05100085	RSO05120085	RSO05150085	RSO05200085
6 mm	RSO06020085	RSO06040085	RSO06060085	RSO06080085	RSO06100085	RSO06120085	RSO06150085	RSO06200085
7 mm	RSO07020085	RSO07040085	RSO07060085	RSO07080085	RSO07100085	RSO07120085	RSO07150085	RSO07200085
8 mm	RSO08020085	RSO08040085	RSO08060085	RSO08080085	RSO08100085	RSO08120085	RSO08150085	RSO08200085
9 mm	–	RSO09040085	RSO09060085	RSO09080085	–	–	–	–
10 mm	–	RSO10040085	RSO10060085	RSO10080085	–	–	–	–
12 mm	–	RSO12040085	RSO12060085	RSO12080085	–	–	–	–

Shaft length: 135 cm

Stent I.D.	Stent length							
	20mm	40mm	60mm	80mm	100mm	120mm	150mm	200mm
5 mm	RSO05020135	RSO05040135	RSO05060135	RSO05080135	RSO05100135	RSO05120135	RSO05150135	RSO05200135
6 mm	RSO06020135	RSO06040135	RSO06060135	RSO06080135	RSO06100135	RSO06120135	RSO06150135	RSO06200135
7 mm	RSO07020135	RSO07040135	RSO07060135	RSO07080135	RSO07100135	RSO07120135	RSO07150135	RSO07200135
8 mm	RSO08020135	RSO08040135	RSO08060135	RSO08080135	RSO08100135	RSO08120135	RSO08150135	RSO08200135

## XI. CE-Certification

CE Marking according	Medical Device Directive 93/42 / EEC, Annex II
Classification of Device	Class IIb
Notified Body	Polish Center for Testing and Certification (PCBC)
CE Number	CE1434
Certificate Number	1434-MDD-301/2020 1434-MDD-301/2020 ANNEX 1
CE-Certificate valid until	2024-05-27