

TFD 7

ALKALINE DETERGENT FOR AUTOMATIC WASHING MACHINES

SURGICAL INSTRUMENTS



AUTOMATIC WASHING MACHINES



■ PRODUCT BENEFITS

- Elimination of different type of secretion, proteins and greases generally dropped on surgical instruments.
- Can be used in ultrasound vats
- Compatibility with a wide range of materials: stainless steel, viton, silicone, PVC, polypropylene...
- Effectiveness guaranteed even in hard water



ALKALINE DETERGENT FOR AUTOMATIC WASHING MACHINES

■ INSTRUCTIONS FOR USE

■ LDI cleaning cycle settings:

- Recommended concentration:.....0.5 to 1%
- Recommended temperature :.....50 to 60°C
- Recommended time contact :.....4 minutes

Concentrations rate can be adapted according to the installation type (washing machine, high-pressure cleaning...) and according to the soil.

Mechanical stirring and heat enhance the detergent action.

■ Product to associate during the cycle:

- Intermediary rinse:
PHOSPHAX (TDS 050.B) :.....acid neutralizer
- Final rinse:
RINCE L7 (TDS 1060.B) :.....rinse additive and neutral active drying
C-NET 2 (TDS 820.B) :.....rinse additive and neutral active drying

■ PRECAUTIONS OF USE

See (SDS 010208)

■ FEATURES

- Form: Liquid
- Color: Colorless to pale yellow
- Pure pH: 12.9 ■ Density: 1,115
- pH at 0.5%: 11.0 ■ pH at 1%: 11.8
- Expiration date: 3 years

■ COMPOSITION

- Sequestrants
- Non-ionic surfactants
- Potassium hydroxide
- Other excipients

■ PACKAGING

- ① Box of 3 low cans x 5 Lref. **1020815**
- ② Box of 4 x 5 L.....ref. **1020820**
- ③ Box of 2 x 10 Lref. **1020822**
- ④ Can of 20 L.....ref. **1020821**



Hazardous product - Follow precautions of use. Always read label and product information before use. Medical Device Class 1.