100

**TearFlo™**

Sterile Tear Flow Test Strips

**TearFlo™**

Sterile Tear Flow Test Strips

CE

2460

100

Ref.: MDT/A1/UU/02 / Issue No.: 01, Dt.: 01/10/2018; Rev No.: 01, Dt.: 01/10/2021

Madhu Instruments Pvt. Ltd.  
A-260, Okhla Industrial Area, Phase-1,  
New Delhi-110 020 India


**Manufactured for and distributed by:**  
MDT J.Zych, A.Budyn sp.j.  
Ul.Skosna 12 A, 30-383 Kraków, Poland  
www.mdt.pl

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### NearFlo™ Sterile Tear Flow Test Strips (Schirmer Test)

**Intended use:** Intended for use during a tear secretion test. The Schirmer test is the oldest method of diagnosis of the dry eye syndrome (DES). **2. Product description** The TearFlow™ strip is a rounded, blunt end that should be placed under the lower eyelid without touching the cornea. The printed scale (0-35 mm), graduated in 2mm increments ( $\pm 0.03$ ), has the accuracy  $\pm 0.1$  mm. Each strip is individually packaged in a sealed and autoclaved medical package and in film to ensure product sterility.

**Indications** As recommended by the ophthalmologist. **4. Contraindications** Do not use in patients with damaged eye surface, inflammation, corneal or conjunctival ulceration, and a case of hypersensitivity to any of the product components. **5. Users** The product is designed exclusively for professional users.

**6. Marking and labeling** The product is marked with the CE mark and having an EU declaration of conformity. **7. Use** Do not apply any topical medication or perform any manipulation on the eye prior to performing the examination. The examination room lighting should be dimmed. Position the patient in the examination chair with the head resting against the headrest. Be careful not to touch the rounded end of the strip when opening the package. With the patient looking up, carefully insert the rounded end under the lower part of the eye lid (about 10mm from the temporal edge or corner of the eye) and then bend the strip down slightly. Wait 5 minutes, remove the strip, and record the length of the moistened area, measuring from the folded end. The strip indicates an initial deficit, and 5-10 mm – an increased deficit, and the aqueous component of tears; <5 mm – advanced DES. **8. Disposal** Because of the possibility of contamination, dispose of with any medical waste. The product is not recyclable.

**R**  
 TearFlow™ Sterile Tear Flow Test Strips – bandelettes de diagnostic  
 (Test de Schirmer)  
 Notice d'utilisation

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