Generic Name: Bilirubin Total Kit (DSA Method)**Abbreviated name: Bil-T (DSA)****Order Information**

Cat. No.	Package size
TBI0102	R1 4×20 mL + R2 1×20 mL
TBI0103	R1 4×32 mL + R2 4×8 mL
TBI0104	R1 4×48 mL + R2 4×12 mL
TBI0105	R1 2×250 mL + R2 1×125 mL

Intended use


In vitro test for the quantitative determination of Bil-T concentration in serum and plasma on photometric systems.

Summary¹

80-85% bilirubin is breakdown product of hemoglobin, and 15-20% bilirubin roots in proteins containing hemoglobin. Measuring the bilirubin in plasma is used for diagnosing and discriminating the reason of jaundice. Hyperbilirubinemia is due to the excess increase of bilirubin and may cause by prehepatic jaundice (e.g. hemolysis), intrahepatic jaundice (e.g. virus hepatitis) or post-hepatic jaundice (e.g. gall-stone). Some chronic and congenital diseases can also result in Hyperbilirubinemia.

Method

Diazotized Sulfanilic Acid (DSA) Method

Reaction Principle

By using a special surfactant to accelerate the solubility of conglutinated bilirubin, total bilirubin with diazo salt at an acid condition to form a red product of azobilirubin. The absorbency increase is directly proportional to the concentration of bilibrubin.

Reagents**Components and concentrations**

R1:	Hydrochloric acid	100 mmol/L
	Sulfanilic acid	5 mmol/L
	Surfactant	1% (m/v)
R2:	Sodium nitrite	72 mmol/L

Warnings and precautions

1. For in vitro diagnostic use only.
2. Take the necessary precautions for the use of laboratory reagents.
3. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
4. Disposal of all waste material should be in accordance with local

guidelines.

5. Material safety data sheet is available for professional user on request.

Reagent Preparation

R1 and R2 are ready to use.

Storage and stability

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

Once opened, the reagents are stable for 40 days at 2-8°C when refrigerated on the analyzer or refrigerator.

Contamination of the reagents must be avoided.

Do not freeze the reagents.

Reagent blank absorbency

The absorbance of reagent blank at 546 nm should be < 0.2 A.

Materials required but not provided

1. Calibrator and controls as indicated below.
2. NaCl solution 9 g/L.
3. General laboratory equipments.

Specimen collection and preparation ²

1. Serum, EDTA plasma is suitable for samples. Whole blood, hemolysis is not recommended for use as a sample. Freshly drawn serum is the preferred specimen.
2. Use the suitable tubes or collection containers and follow the instruction of the manufacturer, avoid effect of the materials of the tubes or other collection containers.
3. Centrifuge samples containing precipitate before performing the assay.
4. Stability:

Store samples protected from light. Perform the assay immediately.

Assay procedure

	Blank	Sample
R1	200 µL	200 µL
Dist. water	3.0 µL	—
Sample	—	3.0 µL
Mix,, incubate for 5 min. at 37°C, read the absorbance, and then add:		
R2	50 µL	50 µL
Mix thoroughly, incubate for 4 min. at 37°C, and read the absorbance again.		
$\Delta A = [\Delta A \text{ sample}] - [\Delta A \text{ blank}]$		

Application sheets for BS series analyzers are available in this document. Refer to the appropriate operator manual for the analyzer-specific assay instructions.

Calibration

1. It is recommended to use the Human multi-calibrator from Mindray and 9 g/L NaCl for two-point calibration. Traceability of the multi-calibrator can refer to the calibrator instructions for use of Mindray Company.

2. Calibration frequency:

After reagent lot changed.

As required following quality control procedures.

Quality control

At least two levels of control material should be analyzed with each batch of samples. In addition, these controls should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or troubleshooting procedures as detailed in the appropriate system manual. We recommend using the Human Assayed Control made by Mindray to verify the performance of the measurement procedure; other suitable control material can be used in addition.

Each laboratory should establish its own internal quality control scheme and procedures for corrective action if controls do not recover within the acceptable tolerances.

Calculation

The analyzer calculates the Bil-T concentration of each sample automatically after calibration.

Conversion factor: $\text{mg/dL} \times 17.1 = \mu\text{mol/L}$

Or: $C_{\text{sample}} = (\Delta A_{\text{sample}} / \Delta A_{\text{calibration}}) \times C_{\text{calibration}}$

Reference Intervals ¹

Each laboratory should establish its own reference intervals based upon its patient population. The reference intervals measured at 37°C listed below were taken from literature:

Sample Type		Conventional Units	S.I. Units
Serum / Plasma	Adult	0.1-1.2 mg/dL	2-21 $\mu\text{mol/L}$

Performance Characteristics

Representative performance data obtained from Mindray system (Mindray BS series analyzers / Mindray Bil-T Reagent) is given below. Results may vary if a different instrument, an individual laboratory or a manual procedure is used.

Limitations-interference

The following substance was tested for interference with this methodology. Criterion: Recovery within $\pm 10\%$ of initial value.

Substance	Level Tested	Observed Effect
Ascorbic acid	10 mg/dL	NSI*
Intralipid	1000 mg/dL	NSI

* NSI: No Significant Interference (within $\pm 10\%$)

Linearity range

The Mindray System provides the following linearity range:

Sample Type	S.I. Units
Serum / Plasma	1.7-600 $\mu\text{mol/L}$

If the value of sample exceeds 600 $\mu\text{mol/L}$, the sample should be diluted with 9 g/L NaCl solution (e.g. 1+ 4) and rerun; the result should be multiplied by 5.

Analytic Sensitivity/Limit of Detection

The lowest measurable Bil-T concentration that can be distinguished from zero is 1.7 $\mu\text{mol/L}$ with 99.7% confidence.

Precision

Precision performance using the CLSI Approved Guideline EP5-A2 to assay serum control appears in the table below³. U: $\mu\text{mol/L}$

Type of	Level I			Level II		
	Mean	SD	CV%	Mean	SD	CV %
Imprecision						
Within-run	18.71	0.22	1.17	92.76	0.94	1.01
Between-run		0.23	1.23		1.33	1.43
Between-day		0.29	1.53		1.49	1.61
Within-device		0.43	2.28		2.21	2.38

Method Comparison

A comparison between Mindray System (Mindray BS series analyzers /Mindray Bil-T Reagent) (y) and Hitachi/Roche System (Hitachi /Roche Bil-T) (x) using 40 samples gave following correlation ($\mu\text{mol/L}$): $y = 1.0226x + 0.0547$, $R^2 = 0.9994$.

Details of the comparison experiments are available on request.

References

- 1.Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 192-202.
- 2.Tietz NW. Clinical Guide to Laboratory Tests.3rd ed. Philadelphia, Pa: WB Saunders Co; 1995:88
- 3.CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

Graphical symbols

In Vitro Diagnostic
medical device

Batch Code

European
Conformity

Authorized representative in
the European Community

Use By

Consult
Instructions for useTemperature
Limit

Manufacturer

Catalogue
number

© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp.
GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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Referências

1. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 192-202.
2. Tietz NW. Clinical Guide to Laboratory Tests. 3rd ed. Philadelphia, Pa: WB Saunders Co; 1995:88
3. CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP5-A2 [ISBN 1-56238-542-9. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087 USA, 2008.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

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© 2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tüm hakları saklıdır
Üretici:Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
AT Temsilcisi:Shanghai International Holding Corp.GmbH(Avrupa)