SIGNA Champion 1.5T

CONFIDENTIAL

SECRET DE AFACERI

DETALII DESPRE PRODUS

I. SIGNA Champion 1.5T 30.1 Winner EMEA

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
CONFIDENTIAL
SECRET DE AFACERI

5941068-1EN
Revision 1
US English

General user documentation.
© 2023 GE HealthCare. GE is a trademark of General Electric Company used under trademark license.
Reproduction and/or distribution is prohibited.

5.7 RF coils

Imaging coils are tuned to match the precessional frequency of nuclei under evaluation. Generally, the length of coil is equal to the FOV the coil covers. The depth of penetration is governed by the coil elements. When selecting a coil, keep in mind the FOV, how deep you need to image, and the size of the patient. Phased array and surface coils need to be placed close to the area of interest. The broad category of imaging coils can be classified into 4 subcategories:

- Head
- Body
- Surface
- Phased array coils


Each coil, other than the body coil, has an operator manual. Refer to the coil operator manual when setting up the patient for an exam.

Head coil

The head coil provides higher SNR (Signal-to-noise ratio) than the **Body** coil due to the smaller size. It is used primarily to image the head, although they can be used for imaging any body part that fits into the coil. It is an example of a volume (uniform depth of signal) coil.

The transmit/receive head coil supplied with your MR system is a Type BF applied part.

Head Coil

CONFIDENTIAL
SECRET DE AFACERI

Body Coil

The **Body Coil** is a transmit/receive volume coil used for large FOV imaging and for uniform depth penetration. The **Body Coil** is located within the magnet enclosure and is invisible to you and the patient. The **Body Coil** can also act as a transmit only coil when used with receive only coils.

CAUTION

Provide all patients with the Patient Alert bulb. This can be especially important for procedures that require the concerted attention of the technologist/operator at the MR or Advantage Workstation (AW) operator console, e.g., fMRI sequences.

The patient alert system is comprised of a squeeze ball that is connected to the PAC and a control box. The control box is typically located at the operator console.

The patient alert port is a Type BF applied part.

Use the following steps to set the patient up with the alert bulb.

Step-by-step instructions

1. Give the patient the alert bulb.

The bulb is a rubber product and not latex.

2. Instruct the patient to loosely hold the alert bulb and to only squeeze the bulb to bring you into the scan room to consult with the patient and attend to his needs.

Figure 6-7 Patient holding alert bulb

CONFIDENTIAL
SECRET DE AFACERI

3. If the patient squeezes the rubber ball end of the alert system, a loud sound is heard in the control room. Follow these steps.
 - 3.1. Press **RESET** to stop the alarm and reactivate it.
 - 3.2. Go into the scan room and attend to the patient's needs.

Figure 6-8 Patient Alert control box. Option 1: Steady/Pulse button on left and Reset button on right. Option 2: Reset button.



Figure 7-17 The invalid state is only confirmed after the series have been saved in RXd state

- If the AutoBind task displays an **INV** (invalid) state and a message displays that the AutoBind requirements are not met, then the series will not bind. The invalid state is only confirmed after the series have been acquired, thus it is very important that all series meet the AutoBind requirements.

10. Click **Save Rx** for each series.
11. Click **Scan** to acquire the data for each series. Alternatively, click **Auto Scan > Scan**.
12. When all the series are completed, and the AutoBind post-process task is in the Done state, select the bound series. From the Workflow Manager control panel, click **View**.

7.1.10 Scan with Auto Paste

Before you begin

- A scan session is open.
- Localizers and other series have been acquired.

CONFIDENTIAL
SECRET DE AFACERI

- For scan details, see [Scan with a Protocol](#) on page 333.

About this task

Use these steps to acquire a scan with a Pasting post-process task. Pasting applications include spine, peripheral vascular and whole body oncology screening exams. The procedure below uses a spine as an example for a pasting exam.


Do not use Pasting post process application with images that demonstrate metal implants.

Step-by-step instructions

1. From the *Protocol* screen, select a protocol for Auto Paste.

Figure 7-18 Example of a whole spine protocol that would use Auto Paste

2. Click the  to send the protocol to the Multi Protocol Basket.
3. From the *Protocol* screen, click **Accept**.
4. From the Workflow Manager, click the protocol for the upper spine range.
5. From the Workflow Manager Control area, click **Add Task > Add Post Processing**.
6. From the *Add Post Processing Task* screen, click **Pasting** for the whole spine application.
7. From the *Add Post Processing Task* screen, click **Setup**.

8. From the *Paste* screen, complete all fields and selections.

Figure 7-19 Paste screen

Paste

Task Description:

Start Processing:

☒ Automatic

☐ Manual

Source Series Data:

Series #	Description
No 1	Type A: 3D Contrast
No 2	Type A: 3D Contrast
No 3	Type A: 3D Contrast

Description:

☒ Equalize Station Intensity

☒ Automatic Registration

☒ Allow Oblique Stations

Arterial

Arterial

Venous

Mask

Accept Cancel

CONFIDENTIAL

SECRET DE AFACERI

- 8.1. Typically select **Automatic**. Alternatively, click **Manual** to paste the images after the scan session has ended
- 8.2. From the Source Series Data list, select the series that you want to Paste. The series should have the same contrast type, for example sagittal T1 FLAIR.
- 8.3. For multi-phase 3D Vascular scans, a menu appears from which you can select Arterial, Venous or Mask to determine the image type for pasting.
- 8.4. From the **Description** text box, type and enter a series description, for example Sagittal T1 spine.
- 8.5. Click **Equalize Station Intensity** option to reduce the difference in the intensity ranges between stations or series and provide more uniform brightness and contrast.
- 8.6. Click **Automatic Registration** option to perform automatic registration of overlap regions.
- 8.7. If needed, click **Allow Oblique Stations** to permit an oblique plane up to 30 degrees between the superior/inferior patient axis.
- 8.8. Click **Accept**
 - The *Paste* screen closes and Auto displays as the status of the post process task in the Workflow Manager.

Figure 7-20 Workflow Manager with Auto task

Task		Series Data		
#	Status	Description		Time
1	Done	1; LOCALIZER		00:16
2	InRx	 Sag T1 FSE UPPER		03:27
2.1	Auto	Auto Paste Sag T1		
3		 Sag T2 FRFSE UPPER	2	03:29
3.1	Auto	Auto Paste Sag T2		
4		Sag T1 FSE LOWER		02:57
5		Sag T2 FRFSE LOWER	4	03:28
6		Ax T2 (ROI)		03:21
7		Ax T1 (ROI)	6	03:02

CONFIDENTIAL
SECRET DE AFACERI

9. To define another Paste task, repeat this step.

10. Graphically prescribe the scan locations.

- For vascular studies, make sure that the vessels are within the same scan plane, to produce optimum Paste images.
- If the Paste task shows INV (Invalid State) and a warning that pasting restrictions are not honored and the images will not be pasted, adjust the scan slices until the Paste task shows Auto. In the Workflow Manager, place your cursor over the INV series to see the roll-over message that describes the reason for the invalid state. The following are several scenarios that result in an Invalid State.
 - If possible, acquire the sagittal and coronal images without an oblique angle.
 - A single oblique axis is allowed. Keep in mind that oblique data sets degrade the accuracy of the pasted series. The automatic and manual Paste tasks display a warning of degraded performance if the image column axis exceeds 0.1 degrees with respect to the superior/inferior axis. The paste task will not complete if the column axis is greater than 30 degrees with respect to the superior/inferior axis. The Paste task will not process series with a dual oblique.

- Always overlap the two series FOV in the scan prescriptions by at least 15% but no greater than 50%. If the overlap is less than 15%, the Registration algorithm does not apply and a prompt appears. If the overlap is greater than 50%, the Paste series is invalid until the overlap is corrected.
11. Click **Save Rx > Scan**.
 12. When the two linked series are completed, and the Paste post-process task is in the Done state, select the pasted series and from the Workflow Manager control panel, click **View**.

7.2 Scan with AutoScan

AutoScan automatically downloads, prescans, and scans all RXD series in the Task List.

AutoScan extends the multi-tasking capabilities by allowing you to prescribe multiple scans in advance. The next RXD series in the list begins automatically when AutoScan is selected. All RXD series in the Workflow Manager list are scanned sequentially.

AutoScan stops if the following conditions exist:

- A table move is required (Move to Scan).
- Move to Scan is not activated after two audio alarms and scan times out.
- A contrast series is in the queue following a non-contrast series. When the first series containing contrast On is ready to scan, the contrast injection message appears. You must click **OK** to proceed. The message does not appear again. AutoScan must be reselected to start AutoScanning mode.
- A new landmark is selected.
- New Patient or a patient from Patient Register list is selected.

Multi-station series can not be AutoScanned.

7.3 Scan with Auto Protocol Optimization

About this task

Use these steps for a quick scan parameter adjustment for breath hold body scans. Multiple optimized scan protocols are offered, which have different scan time, resolution.

NOTE

Auto Protocol Optimization supports abdominal Anatomical Regions.

Step-by-step instructions

1. Open a scan session and load one of the following body breath hold series as an InRx state:
 - 3D applications: LAVA, LAVA-Flex, Dual Echo, FRFSE-XL with MRCP Imaging Option, IDEAL IQ, DISCO
 - 2D applications: Fast SPGR, Dual Echo, FIESTA, SSFSE, FRFSE-XL, FSE-XL, DWI EPI
2. From the *Scan Parameters* screen, click **GRx** and set up the scan locations on the localizer images.
3. From the , click .
 - The button is activated once all parameters are filled and the **Save Rx** button is activated.

- The Image Management work area. From the patient list, press **Ctrl** and simultaneously click the pre-contrast series and the post-contrast series to select them both and then click **READY View** from the Session Management list.
- From the READY View protocol screen, click **Ser** and follow the instructions detailed in the READY View section.
- VIBRANT images can be post-processed using the purchasable CADstream option. For details, see the CADstream operator manual.

 WARNING

It is possible that a spatial distortion can be seen on 3D data sets, especially in the lateral-most VIBRANT images. The distortion can be demonstrated in sagittal versus axial data sets. There is a potential risk for lesion localization misregistration during biopsy procedures, which could result in a re-biopsy of the patient.

 CAUTION

Images labeled as water may include signal from fatty tissue, and images labeled as fat may include signal from water. This error may occur in regions of high magnetic field variation, in spatially isolated tissue, due to patient or tissue motion, due to phase wrap artifacts, due to similar water and fat intensity in T2-weighted Flex scans, due to TE values beyond recommended limits, and/or in images with low signal-to-noise ratios. The presence of fat tissue in images labeled as water, or vice versa, may occur within single images or throughout an entire stack of slices. By default, both sets of images (labeled fat and labeled water) will be reconstructed and inserted into the database for review. Proper calibration and center frequency selection will reduce the occurrence of this error. Complete elimination of this error may not be possible and thus interpretation of MR images must be completed by trained personnel.

Results

Image annotation

M3D/VIBRANT/flip angle

Water and Fat images are synthesized from collected in-phase and out-of-phase images; thus TE values for Water and Fat images are the average of IN-phase and Out-of-Phase TEs. In-phase and out-of-phase gets annotated according to the TE and TE2 values on screen.

8.2 Scan Anatomy: Cardiac

8.2.1 Cardiac introduction

Cardiac Gating/Triggering allows the same data to be collected at the same points of the cardiac cycle for each repetition of the pulse sequence allowing motion artifacts to be reduced, or sometimes eliminated. The benefit is that the motion is always the same and is essentially “frozen,” allowing motion artifacts to be reduced or, sometimes, eliminated. Gating and triggering of MR images can be performed with ECG or Peripheral gating techniques.

8.2.1.9 ECG waveform

Electrical impulses cause the heart to contract and, therefore, blood to flow throughout the body. The electrical activity of the heart can be detected by measuring the voltage difference between electrodes attached to the patient. The voltage differences can be mapped by an ECG and the resulting ECG waveform can be used during cardiac gating/triggering to reduce pulsatile cardiac motion by synchronizing data acquisition to the cardiac cycle. The ECG gating devices are safe for patients due to the fiber-optic technology.

ECG gating can eliminate motion artifacts by using the heart's electrical activity as reported by the ECG to trigger the data acquisition. It is useful during the following applications:

- Cine-PC imaging to examine flow
- Heart imaging for structure or function
- Thorax imaging

Use these steps to prepare your patient for an ECG gated cardiac examination. To perform a successful gated examination, it is very important to follow these steps for proper ECG lead placement and cable routing. Attach electrodes no more than 5 to 10 minutes before a scan.

The following are Type BF applied parts:

- ECG lead connector
- respiratory bellows connector

The ECG waveform represents the electrical activity of the heart that correlates to the following heart motion: cardiac contraction called systole, and cardiac relaxation called diastole. The figure below displays a normal ECG pattern of the waves in the cardiac cycle.

CONFIDENTIAL
SECRET DE AFACERI

Graphic Rx

Graphic Rx is available whenever a valid localizer image exists. You can graphically prescribe an ROI with the MRS sequences. Graphic ROI prescription requires that the scan plane of localizer image(s) must be the same as the prescribed scan plane. Click on the image to display the box cursor (ROI).

Phase cycling

Phase cycling is the repetition of a pulse sequence and signal acquisition such that all acquisition parameters except the phase of the RF pulses (and, possibly, the receiver phase) remain unchanged from repetition to repetition. Phase cycling suppresses or eliminates undesirable signals while taking advantage of the effects of signal averaging. In the current spectroscopy pulse sequences, only the RF phase is changed, or cycled, in the repetitions.

In one of the simpler examples of a phase cycling scheme, a pair of excitations is obtained with the phase of an RF pulse used in the second acquisition inverted relative to the phase of the RF pulse used in the first acquisition, causing an inversion of the desired signal between the two acquisitions. The second excitation is then subtracted from the first to produce a single frame. The result is that the first signal and the inverted second signal are added together while the RF errors independent of the phase of the RF pulses are canceled by the subtraction.

The spectroscopy sequences always use the largest phase cycle that is compatible with the selected NEX value.

Homogeneity adjustments - gradient shimming

It is possible to manually adjust the linear x, y, and z shims by adjusting the corresponding gradient current offsets on the Spectroscopy screen. However, adjusting the homogeneity through manual shimming is neither fast nor easy, particularly for nuclei other than hydrogen. For most users and for most applications a better choice is to rely on the Autoslim capabilities of the system.

The hydrogen based Autoslim capability is both fast and reliable, and should be used to optimize the shim whenever possible.

Generally, Autoslim should be selected during the acquisition of the localizer image(s).

For voxel localized spectroscopy, the shim through the voxel can be improved by running Auto Prescan with Autoslim for the slice containing the voxel prior to prescribing the spectroscopy data acquisition.

For slice localized spectroscopy, the shim across the slice can be optimized by running Autoslim as part of an imaging sequence Auto Prescan process for the same slice that is selected in the spectroscopy data prescription. If the focus of the acquisition is a small region of the slice, prescribe a PRESS CSI scan with a voxel that covers the small region, and use Auto Prescan with Autoslim to optimize the shim through the voxel.

Spectroscopy data

The basic units of spectroscopy data are the complex point, the excitation, and the frame.

- **Complex Point:** The complex point is the simplest data unit. It is produced from the signals acquired by the receiver during data acquisition. The digitized signals are stored as 4 byte integers, or, if the Extended Dynamic Range Imaging Option is selected, 8 byte integers. According to standard imaging usage, the complex point consists of an I and Q pair, where I is the in-phase component, and Q is the quadrature component. The I and Q components are treated as the real and imaginary parts, respectively, of a complex number.
- **Excitation:** An excitation is the collection of a selected number of complex points acquired during a single pass through a pulse sequence. For spectroscopic data acquisitions, an excitation consists

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AW Family

● AW VolumeShare 7

The AW Family is GE Healthcare's ecosystem of processing and advanced visualization solutions across platforms and applications.

Advantage Workstation (AW) with VolumeShare 7 enables multi-modality image review, comparison, and processing with simplicity and power. Featuring, 64-bit technology, this workstation is a cornerstone solution to Radiology departments around the world.

This document is intended to overview the features, specifications, use cases, and other key information of AW VolumeShare 7.

CONFIDENTIAL
SECRET DE AFACERI

Contents

Overview.....3

 Key Features 3

 Industry Standards 3

 DICOM Conformance Standards: 3

 Filming Protocols 3

 Indications for Use 4

 Regulatory Compliance 4

Product Details5

 Modalities 5

 License Management 5

 User Interface 5

 Standard Applications & Capabilities 6

 Workflow Management 9

 Security Capabilities 11

System Components 11

 Workstation Configuration 11

About GE HealthCare 13

CONFIDENTIAL

SECRET DE AFACERI

Overview

Key Features

- Simple drag and drop action for networking, media interchange, and filming.
- Multi-tasking capabilities.
- Active directory integration allows enterprise level user authentication. Directories supported include Microsoft® Active Directory® (MSAD) and other LDAP authentication services.
- Patient list to manage images from local workstation or remote systems.
- Postfetch feature, designed to optimize your reading workflow by automatically gathering a patient's prior exams according to your pre-set criteria.
- Support for Key Images.
- End Review automates routine filming and networking tasks with just one click.
- Search Advantage for fast and easy search of a patient's exam history on PACS or any other DICOM®-compliant device.
- Enhanced Quick Filters of the Patient List filters studies by Modality, Date, End Review status or Exam Description.
- DICOM CD/DVD/USB Creation Tool.
- 2D Viewer for image display, manipulation, annotation, review.
- Integrated Filmer with enhanced flexibility to perform all filming and data exporting tasks.
- Access to a wide variety of applications for greater diagnostic flexibility.

Industry Standards

AW VolumeShare 7 complies with a wide variety of industry standards to facilitate adoption of features and performance improvements as the computing and medical imaging industry evolves.

DICOM Conformance Standards:

- DICOM 3.0 Storage Service Class for RT, CT, MR, CR, X-ray (Angio and R&F), Digital X-ray (DX), MG, NM, PET, U/S, Secondary Capture, Secondary Capture Color DICOM Image Objects. (Service Class User (SCU) for image send and Service Class Provider (SCP) for image receive).
- DICOM 3.0 Query/Retrieve Service Class (SCU and SCP).
- DICOM 3.0 Storage Commitment Service Class (SCU).
- DICOM Print (Color and B&W).
- DICOM Media Interchange (CD-R, DVD+R(W)).

Filming Protocols

DICOM Print (Color and B&W) and Adobe® Postscript (Color and B&W) for supported printers.

CONFIDENTIAL
SECRET DE AFACERI

Indications for Use

AW VolumeShare 7 is a review workstation, which allows easy selection, review, processing and filming of multi-modality DICOM images from a variety of diagnostic imaging systems. When interpreted by a trained physician, filmed or displayed images on the AW monitor may be used as a basis for diagnosis, except in the case of mammography images.

Regulatory Compliance

This product complies with the Regulation (EU) 2017/745 of the European Parliament and the Council on the medical device (MDR).

This product or its features may not be available in some countries or regions. Please contact your sales associate.

Rx Only

CONFIDENTIAL
SECRET DE AFACERI

Product Details

The Advantage Workstation VolumeShare 7 software technology makes a multi-modality advanced visualization workflow solution that helps to enhance diagnostic confidence and productivity. This solution features software that is optimized for 64-bit technology and multicore processor hardware to provide leading edge performance. It includes AW's premier 3D image analysis package, Volume Viewer, with a suite of volumetric visualization, and analysis tools for CT, MR, 3D X-ray and PET.

In addition, this package includes the multi-modality 2D Viewer with dedicated functions for review of CT, MR, X-Ray, Angio, DX, U/S and PET images. It also includes the Filmer, a multimedia export tool for creating electronic films, filming, and exporting to an internal web server, CD, or DVD.

AW VolumeShare 7 is also available for purchase without Volume Viewer for users who require only 2D image applications. Volume Viewer dependent features will not be supported with this configuration.

Modalities

Supported DICOM modalities include:

- Computed tomography (CT),
- Magnetic resonance imaging (MR),
- Radiofluoroscopy (RF),
- Xray angiography (XA),
- Computed radiography (CR),
- Digital radiography (DX),
- Mammography (MG),

- Nuclear medicine (NM),
- Positron emission tomography (PET),
- Ultrasound (US)
- Structured Reporting (SR),
- Key Objects (KO)

License Management

- Application Usage Monitor tracks how often a particular license was used, when, and by whom. You can compile and view usage reports for a specific period of time. This helps you understand how your systems contribute to departmental productivity and to plan for future license requirements based on actual usage data.
- For floating licenses, a secondary license server can be assigned for redundancy. If the primary license server becomes unavailable, the secondary can be contacted so that work can continue.

User Interface

- Displays on one or two 1280x1024 monitors with scroll or optical mouse.
- Simple drag-and-drop mechanism for networking, media interchange, filming and data export options.
- Automatic help message display briefly describes function currently under the mouse pointer.
- Shortcut keys and programmable window/level function keys for accelerated control.
- Multi-tasking capabilities that allow Patient List, 2D Viewer, Filmer and Volume Viewer to run simultaneously with Fast Switch capabilities between applications.

CONFIDENTIAL
SECRET DE AFACERI

Standard Applications & Capabilities

Separate application licenses must be purchased to access advanced applications functionality.

2D Viewer

The 2D Viewer is an application used to display, manipulate, annotate and review 2D images by a trained physician for diagnostic interpretations.

Display Customization

Display customization allows the user to manage layouts to display data. Key benefits include:

- Allows modification of number of exams or series displayed to facilitate either single or multi-exam reviews.
- Standard layouts provide flexibility to tailor the image display from 1x1 to 8x8. The current layout is retained if the images are from the same modality and orientation.
- Annotation levels allow selection of the image information fields to display.
- Toolbar customization enables control of which buttons are displayed based on user preference.

Study Navigation

- The navigator lets users assign a series to a view on the fly.
- Cine mode also provides temporal, spatial, or manual playback loops.
- With two exams or series loaded, Cine mode supports a side-by-side display format with synchronized playback loops for more efficient comparisons.

Image Review

The initial image window and level setting is based on the DICOM header. Once displayed, several methods of adjusting image window and level are provided.

Provides routine image manipulation features:

- Flip/Rotate, Zoom, Pan, Magnifying Glass,
- Inverse Video: inverts grayscale color map,
- Display normal: lets you return the image to its default viewing parameters.

You can access the following features using a single mouse click directly on the image. Having these direct manipulation tools easily at hand gives you fast access with less distraction from your review task by eliminating the need to return to the graphical user interface controls: 2D distance, Angle, Report cursor, Box ROI, elliptical ROI and free-hand ROI.

The program continuously updates statistics on the fly.

- The Annotation feature lets you highlight areas of interest by adding text and line/arrow pointers to any image feature.
- The Copy/Paste/Erase feature lets you copy/paste/remove any text or graphic placed on an image.
- Cross-reference indicates the position of the current slice over the scout or localizer or any other non-parallel series.
- Save lets you store a copy of the image as it appears on the screen for future review.
- The Key Images feature allows you to flag images and create key objects.

Advanced X-Ray Analysis

The user can apply shutters to X-ray images to focus on specific areas within the image. Shutters are black, opaque overlays with elliptical or

CONFIDENTIAL
SECRET DE AFACERI

rectangular cutouts that you can size as desired. Once a shutter is applied, you can move the underlying image to shift the viewable area. Shutters can also be applied to CT and MR images with the Image Matte feature.

The user can apply several levels of edge enhancement filters to an image.

- **Image Subtraction:** A graphical user interface lets you select a mask and subtraction of all images in an associated sequence. Subtracted images may be saved as a sequence with one mouse click.
- **Landscape:** With this tool, you can introduce a percentage of the mask image into the subtracted image for anatomical reference. A graphical user interface lets you specify the percentage of mask to be reintroduced.
- **Pixel Shift:** Auto pixel shift optimizes a mask to image registration. You can also move a mask manually to optimize subtraction quality in a specific region of interest.
- **Split Pixel Shift:** You can split the screen horizontally or vertically for pixel shift.
- **Maximum/Minimum Opacifications:** You can integrate selected images to provide a resulting Max. Op./Min. Op. image.

Batch Filming

The Print Series feature lets users automatically batch film an entire series with a single keystroke.

Volume Viewer

Volume Viewer is a multi-modality AW Family advanced visualization application. It provides excellent 3D visualization and processing capabilities for reading and comparing CT, MR, 3D X-Ray, PET, PET/MR, and PET/CT datasets. It

features a breadth of high-performance analysis tools such as Multiplanar Reformat, Curvilinear Reformatting, Volume Rendering and Lumen Navigation.

Volume Viewer is the pre-requisite application for what are known as Volume Viewer Foundation (VVF) applications within the AW Family Advanced Visualization Application portfolio. VVF applications share a similar user interface and have some interoperability. For detailed description of Volume Viewer features, please refer to the Volume Viewer Product Data Sheet.

Filmer

The integrated Filmer enhances the efficiency of the review station and gives users greater filming and exporting flexibility. This feature supports two modes: the Mini Filmer mode and Full-Screen mode, which provides the ability for customization and film layout template creation. Three key mechanisms give users flexibility:

- Free format filming
- DICOM structured reporting (SR)
- Data export (HTML/PDF and JPEG, PNG, MPEG, AVI, or QTVR)

With the Filmer, users can easily extract significant images from any AW Family AV application (2D Viewer, Volume Viewer, Ready View, etc.).

- You can film images individually by dragging and dropping to the on-screen Filmer, or by a single F1 keystroke.
- Multiple image formatting allows filming multiple images in a single page frame with the F2 keystroke.
- With Film MID users can send multiple images to a single Filmer frame.

CONFIDENTIAL
SECRET DE AFACERI

- Batch filming is supported by applications that provide that capability (Print Series in the 2D Viewer, Batch Film Protocols in Volume Viewer).

Communication between applications and the Filmer is accomplished in the Mini Filmer mode, which provides a minimized footprint. The Mini Filmer mode provides the following:

- Compatibility with Batch Filming from Volume Viewer (may be optional in your configuration)
- Store and position images transferred from an application.
- Type of export (film, media, database).
- Rapid switching between Full-Screen Filmer and application.

The resultant electronic films become a quick summary of the patient study and radiological interpretation which can then be reviewed by clinicians and physicians. They can also be saved as independent files (DICOM SR and Secondary Captures) for teaching purposes. Electronic film can contain one or several pages with specific layouts for each page or all pages.

A flexible Edit Mode provides the ability to easily add, manipulate, format, or delete images from the film. Images can contain text and graphics from measurements and your annotations, and may be window/leveled, magnified, flipped, rotated, or cine. Additional annotations can be added to the image in edit mode.

Preview Mode displays the film as it will be printed or exported, considering the layout applied to each page and the compression level specified for non-DICOM exporting (JPEG/PNG and MPEG/AVI/QTVR).

Printing

AW VolumeShare 7 includes Network DICOM Print (B&W and color). For supported printers, AW VolumeShare 7 includes Network PostScript capability (B&W and color). For PostScript and DICOM printers, any printing format created in the Filmer is supported (e.g., non-square matrix formats for rectangular images such as CT run-offs).

- AW VolumeShare 7 postscript printing has been tested on the following devices: Codonics 1660M, 1660MD or Horizon, Lexmark Optra 1650N, 1855N, SC1275N, C710N, C720N, T612 or T614, Seiko 1720D, Kodak DMI3600, Quantum GL2101HD, declared with film/thick paper or with plain paper, Tally T8106, HP LaserJet and Xerox Phaser
- Digital cameras and analog cameras using a 3M-952 protocol (including DASM interface) are not supported on AW VolumeShare 7.

Exporting

Data export is integrated in the Filmer, thus providing all image processing tools needed for multimedia image export. **The Filmer exports any electronic film to DICOM SR, PDF/HTML, JPEG, PNG, MPEG, AVI or QTVR format.** AW VolumeShare 7 enhanced security of the export functionality with the use of SFTP connectivity. Non-DICOM data can also be saved on a multi-session CD/DVD. The Data Export capability is intended only for publishing and communication, not for diagnostic purposes. Its simplicity is reflected in the different export mechanisms available:

- CD/DVD removable media
- USB flash drive
- Network HTTP and FTP protocols

CONFIDENTIAL
SECRET DE AFACERI

Cardiac Review and Export

Processing and review of CT, MR and PET cardiac exams with manual oblique reformatted protocols can be exported as a multi-phase Cine movie that allows referring physicians to review exams in a dynamic mode.

Workflow Management

AW VolumeShare 7 provides features designed to optimize your reading workflow:

- The Patient List provides tools that enable, you to sort and filter imaging studies stored on the local workstation or on remote systems.
- Postfetch allows you to retrieve a patient's prior related DICOM exams from a remote DICOM host. Retrieval can be triggered by the arrival of a new patient study on the workstation from the network or supported media.
- Support for IHE Key Image Notes profile, allowing you to mark certain images in the 2D or 3D viewer as key images. Key images are displayed as a separate series in the exam list and accessed or sent to an IHE KIN-supporting PACS. Specific tags and Key Image Notes may be created using the 2D Viewer application.
- Support for external USB disks to serve as DICOM storage media. The USB disk can be accessed for reading and writing via the media button, the same way you would access a CD or DVD. The maximum number of images that can be stored on a USB device is restricted to 300,000 across all exams. USB media must be formatted using the FAT32 file system. See the AW VolumeShare 7 Basic Display/Viewer/Filmer User Manual for more details on USB storage.

- Preference sharing among users. Sharable preferences include: Filmer layouts, Volume Viewer custom protocols, and viewer preferences. A standard USB flash drive may also be used to share preferences between workstations.
- One-Touch protocols enable you to define an application or preset protocol to launch automatically based on DICOM elements.
- DICOM Query/Retrieve Storage Class User (SCU) and Storage Class Provider (SCP) provide seamless network integration.
- DICOM Storage Commitment SCU lets you know when exams have been archived successfully on DICOM devices, such as PACS, that support Storage Commitment SCP.

End Review

End Review automates the routine tasks required at the end of reviewing each exam. The “End Review” flag in the patient list allows you to mark exams as “Done” after post-processing has been completed. End Review allows you to automatically perform one or several of these actions:

- Print pages prepared in the Filmer to the default printer (DICOM or Postscript filmers supported), with the option to automatically clear the Filmer after printing.
- Save electronic films in the workstation's attached DICOM database.
- Automatically push the entire exam or only the series created on the AW to one or more remote hosts.

Remote Network Host Management

AW VolumeShare 7 Patient List provides a control panel to select a remote host or destination for

CONFIDENTIAL
SECRET DE AFACERI

network transfer. Icons signify different device types (acquisition, post processing, PACS, etc.) that are accessible to the workstation via the network.

You can initiate a network transfer simply by dragging an exam, patient, series, or image(s) to the remote host icon. You can activate a remote browser to display detailed exam information from the remote host, including filtering information if supported.

Search Advantage

This advanced search engine lets you search a patient's exam history on PACS or other DICOM device with a few mouse clicks.

Worklists and Filters

You can filter the patient list by modality, date, end review status, or exam description. Most filters are available on the remote host patient list as well.

You can filter the patient list further by choosing one or any of these parameters:

- Modality.
- Patient name or patient ID.
- Exam location (hospital name) and exam description.
- Series description
- Date and time of day: today, or specified date or date ranges with a specified time or time ranges.
- Radiologist's or referring physician's name.
- Accession number.

Quick Sort (ascending or descending) and Quick Access (entry field) are available for any of the fields displayed at exam level (patient name, patient ID, exam location and description,

radiologist's and referring physician's name, date and time).

Queues

AW VolumeShare 7 manages three queues: network transfer; media interchange, filming. A menu on the patient list gives you easy access to queue status, and an animated icon gives you continuous network queue status information.

DICOM CD/DVD/USB Creation Tool

- The DICOM media creation tool offers you the flexibility to save and retrieve from supported CD, DVD, or USB storage devices.
- Add selected exam, series, or images by dragging and dropping to the pre-mastering window.
- CD and DVD usage percent is shown before you launch the writing process.
- You can optionally down sample certain 1024²x512² X-ray angiographic images during CD creation.
- CD/DVD composer lets you manage large amounts of data. If data size is larger than medium capacity, you are prompted to supply additional media storage space.
- You can specify the number of copies created during a CD/DVD save session.
- You can embed a DICOM viewer on DICOM CD/DVD so data can be reviewed on a PC running a Windows operating system.
- You get lossless JPEG reading of any DICOM CD/DVD media. Lossless writing is available for X-ray exams.
- CD/DVD drives operate at 16x write speeds and read speeds up to 48x depending on choice of media.

CONFIDENTIAL

SECRET DE AFACERI

- For limitations on USB storage, refer to the AW VolumeShare 7 Basic Display/Viewer/Filmer User Manual.

Database Management

- The database management system classifies data according to patient folder description of the DICOM standard: Patient/Exam; Study/Series; Sequence/Images.
- Automatically deletes images on first-in first-out (FIFO) basis, which can be toggled on/off.
- Lock Exam feature protects specified exams from deletion.
- Permanent display of available disk space facilitates disk management.
- Study anonymization tool modifies DICOM elements to remove protected health information.

Security Capabilities

- Full support of Enterprise Authentication with Microsoft® Active Directory® and Novell® eDirectory™. Users can access their AW Workstation with existing Single Sign On credentials. This may streamline HIPAA compliance.
- If Enterprise Authentication is not used, an administrator-defined UNIX login and password is required per workstation. A lock screen feature enables you to lock the workstation display when leaving it momentarily, preventing unauthorized access to patient data.
- Main actions performed on the AW (save, film, network, export) are flagged and saved to an audit trail which includes username, ID, exam ID, date, and time.
- Multiple users can have different local UNIX logins, use the same AW workstation, and yet keep their user preferences unique.

- AW VolumeShare 7 is secured by Product Network Filters (PNF), which keeps unauthorized users from accessing open ports over the network. You can configure this feature to provide only specific remote devices access to AW ports.
- While not recommended, AW VolumeShare 7 does support ClamAV® anti-virus protection. This program is not activated by default in order to guarantee optimal performance of the AW workstation. If needed, your GE Healthcare field engineer can activate the program. However, it is the customer's responsibility to provide a safe internet connection from the workstation and to keep virus definitions current.
- AW VolumeShare 7 provides de-identification and encryption capabilities to limit privacy risks to sensitive information. The patient data exported during clinical workflow may be encrypted in flight by using either the SFTP or DICOM TLS protocol depending on the required dataflow.

System Components

Workstation Configuration

- HP Z4 G4 Workstation
- SLES 15 operating system
- Intel® Xeon® W-2245 Eight physical core/ sixteen logical cores 3.9 GHz CPU
- 64GB (4x16GB) DDR4 2933 MHz or higher ECC Registered DIMM
- NVIDIA Quadro NVS P620, 2GB Graphics card
- 1 x 512GB M.2 NVMe Drive for OS and Applications
- 1 x 1TB M.2 NVMe Drive for image cache.
- Cache, which is subject to overhead can store approximately:
 - 6,400,000 256*256*16 bits uncompressed images OR

CONFIDENTIAL

SECRET DE AFACERI

- 1,600,000 512*512*16 bits uncompressed images OR
- 400,000 1024*1024*16 bits uncompressed images OR
- 100,000 2048*2048*16 bits uncompressed images
- Internal DVD Writer drive for read/write of DICOM CD/DVD media, read/write of Data Export CD/DVD data and service use (DVD Install)
- Integrated dual Ethernet 10/100/1000 Mbit/s ports.
- Additional single Ethernet 10/100/1000 Mbit/s port
- 1 USB QWERTY (or regional) Keyboard and Mouse

Footprint

- Height 38.60 cm (15.2 in.)
- Width 16.89 cm (6.65 in.)
- Depth 44.47 cm (17.5 in.)
- Approximate Weight 17 kg (38 lbs.)

Operating Environment

- Temperature: +5°C to +35°C (40° to 95° F)
- Humidity: 10% to 85% (relative non-condensing)
- Altitude: 0 to 5000m (16,404 ft.)
- Acoustics: LWAd less than 4.3 Bels
- Shock: 40 G peak, half-sine, 2-3 ms

Non-Operating Environment

- Temperature: -40°C to +60°C (-40° to 140° F)
- Humidity: 10% to 90% (relative, non-condensing)
- Altitude: 0 to 12192m (37,000 ft.)

Monitors

- (2) 19" color Flat Panel LCD monitors

- DICOM Part 14 factory calibrated
- Native resolution 1280 x 1024 (5:4 aspect ratio)
- Approximate Weight: 5.3 Kg (11.68 lbs.).
- AC 100 - 240 V : 50 / 60 Hz

Image Networking

- Standard 10/100/1000 Base-T Ethernet for DICOM
- 1000 Base-T dedicated network for optimal Direct Connect performance
- Protocols supported:
 - DICOM 3.0 Storage SCU/SCP and Query/Retrieve SCU/SCP
 - RSVP
 - TCP/IP network layer
- SdCNet is no longer supported
- AW VolumeShare 7 does not support the AdvantageNET network protocol.
- AW VolumeShare 7 does not support DICOM images from GE Healthcare Signa™ version 5.x 1.5T MR systems
- AW VolumeShare 7 software is supported on the previous generation AW VolumeShare 5 HP Z820 and Z440 workstations.

CONFIDENTIAL
SECRET DE AFACERI

GE HealthCare

About GE HealthCare

GE HealthCare is a leading global medical technology, pharmaceutical diagnostics, and digital solutions innovator, dedicated to providing integrated solutions, services, and data analytics to make hospitals more efficient, clinicians more effective, therapies more precise, and patients healthier and happier. Serving patients and providers for more than 100 years, GE HealthCare is advancing personalized, connected, and compassionate care, while simplifying the patient's journey across the care pathway. Together our Imaging, Ultrasound, Patient Care Solutions, and Pharmaceutical Diagnostics businesses help improve patient care from prevention and screening, to diagnosis, treatment, therapy, and monitoring. We are an \$18 billion business with 51,000 employees working to create a world where healthcare has no limits.

Follow us on [Facebook](#), [LinkedIn](#), [Twitter](#), [Instagram](#) and [Insights](#) for the latest news, or visit our website [gehealthcare.com](#) for more information.

Products mentioned in the material may be subject to government regulations and may not be available in all countries. Shipment and effective sale can only occur after approval from the regulator. Please check with local GE HealthCare representative for details.

GE and the GE monogram are trademarks of General Electric Company. Used under trademark license.

Adobe is either a registered trademark or trademark of Adobe Systems Incorporated in the United States and/or other countries.

ClamAV is a registered trademark of Cisco Systems.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

Intel and Xeon, are either registered trademarks or trademarks of Intel Corporation in the United States and/or other countries.

Microsoft, Active Directory, and Windows, are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Novell and eDirectory are registered trademarks of Novell, Inc. in the United States and other countries.

Oracle and Java are registered trademarks of Oracle and/or its affiliates.

All other product names and logos are trademarks or registered trademarks of their respective companies.

©2023 GE HealthCare.

September 2023
DOC1845651

CONFIDENTIAL
SECRET DE AFACERI