Cl⁻ Electrode

Chloride Electrode

(For use with direct ISE module)

INSTRUCTIONS FOR USE

mindray

Cl⁻ Electrode**Chloride Electrode****(For use with direct ISE module)****Order Information**

Cat. No.	Package size
115-084091-00	1 electrode

Intended Purpose

Chloride electrode is used for the in vitro quantitative determination of Cl⁻ concentration in human serum, plasma and urine on Mindray BS series chemistry analyzers. It is intended to be used for diagnosing and monitoring the therapeutic effect of hyperchloremia and hypochloremia and aiding to diagnose renal diseases.

Summary

It is used to quantitatively measure the concentration of Cl⁻ in serum, plasma and urine. Cl⁻ is the main anion of extracellular fluid. Its physiological function is to maintain the metabolism of cells, adjust the osmotic pressure and acid-base balance, playing an important role in maintaining the normal metabolism of nerve and muscle cells. Almost all the chloride ions are absorbed from the intestinal tract and excreted through the kidneys. The detection results of Cl⁻ are used to diagnose hypochloremia and hyperchloremia. Hypochloremia usually occurs in patients with decreased intake, anorexia, coma, digestive tract disease, inability to eat after major gastrointestinal surgery, or adrenal cortex disease and renal disease.

Hyperchloremia often occurs in dehydration, acute nephritis, hypertension, renal failure, or high dietary or parenteral chloride intake.¹

Assay Principle

The sample concentration of Cl^- in serum, plasma or urine on Mindray chemistry analyzer is measured by ion-selective electrode method (direct method). During the measurement of $\text{Na}^+/\text{K}^+/\text{Cl}^-$, the sample flows through the ion measuring electrode and the reference electrode. Due to characteristics of the electrode membrane, the ion to be measured will undergo ion exchange and diffusion on the electrode membrane material to form membrane potential. The response intensity of membrane potential is related to the concentration of the ion to be measured. The higher the concentration of the ion to be measured, the greater the membrane potential. The relationship between the potential of the Cl^- ion selective electrode and the concentration of the Cl^- ion in the sample follows Nernst equation. The concentration of the ion in the sample can be calculated based on the measured potential.²

Key Components

The main components of Cl^- electrode are the sensitive membrane in contact with the sample flow channel and the reference solution and reference electrode Ag/AgCl in the carrier. The components of Cl^- electrode sensitive membrane are PVCs.

Storage and stability

Storage and stability: Refer to expiration date indicated on the label

when stored unopened for 12 months at 2-32°C and without condensation.

On-board stability: After the electrodes are installed, 10,000 samples can be tested for six months. The expiration date is determined by the date when the test period or tested sample volume reaches the limit first.

For replacement, refer to instructions on the Mindray Chemistry Analyzer Operator's Manual.

Specimen collection and preparation

■ Specimen types

Serum, plasma (only heparin lithium anticoagulant is accepted) and urine (to collect urine without additives within 24 hours, remove turbidity from urine by centrifugation, and do not acidify.) are suitable for samples.³

Blood collection tubes from different manufacturers may contain different substances, which may affect the test results in some cases. Each laboratory should establish its own criteria to determine the applicability of collection tubes.

■ Preparation for Analysis

1. Specimens must be centrifuged according to tube manufacturer's specifications.
2. Specimens should be tested as soon as possible after sample collection and pre-analytical treatment.
3. Serum or plasma samples used for measuring Cl^- concentration must be collected in a proper way to minimize hemolysis.

Contamination, clots, and floccules must be avoided.

4. Detection results may be affected if there are such interferents as hemolysis, lipemia and fibrin in the sample. In the above situations, re-collection is recommended. Do not use contaminated specimens.

5. Before loading the specimens on the analyzer, make sure there are no air bubbles.

■ Specimen Storage Stability

The stability of Cl⁻ in serum, plasma and urine specimens stored in sealed test tubes (the deviation index of specimen stability is $\pm 10\%$) is as follows:

Cl⁻ in serum specimens: 7 days at room temperature (15-25°C), 7 days at refrigeration temperature (2-8°C), and 1 year at (-25)-(-15)°C.^{4,5,11}

Cl⁻ in urine specimens: 7 days at room temperature (15-25°C), 7 days refrigeration temperature (2-8°C), and 7 days at (-25)-(-15)°C.^{4,5,11}

Samples cannot be frozen and thawed repeatedly during storage.

For longer storage periods, each laboratory should establish its own sample stability criteria.

Preparation

ISE Reagent Pack and K⁺/Na⁺/Ref electrode are ready to use.

Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

Materials required but not provided

Cat. No.	Product Description	Package size
----------	---------------------	--------------

115-084089-00	Potassium Electrode	1 electrode
115-084090-00	Sodium Electrode	1 electrode
115-084088-00	Reference Electrode	1 electrode
105-025539-00	ISE Reagent Pack	Buffer A: 460 mL Buffer B: 360 mL
105-025540-00	Urine Diluent	1x50 mL
105-025543-00	ISE Cleaning Solution	1x50 mL
105-041321-A0	PROBE CLEANSER	1x50 mL
105-002225-A0	PROBE CLEANSER	1x50 mL
105-025028-A0	DETERGENT C	6x20 mL

For the PROBE CLEANSER and DETERGENT C, only the names and specifications are different. The formula and use are the same.

Calibrator and Control: Please check the section of reagent instruction of Calibration and Quality Control.

Please refer to the manual for information about storage and stability.

For information on materials required for operation of the instrument, and maintenance procedures refer to the Mindray Chemistry Analyzer Operator's Manual.

Assay procedure

Parameters Item		BS-600M chemistry analyzer
Sample types	Serum, plasma	Urine: Dilution first. Mix sample with urine diluent with a proportion of 1:9.

Parameters Item**BS-600M chemistry analyzer**

Sample	Serum, plasma	70 μ L
	Urine	140 μ L

For Mindray BS series chemistry analyzers, the Parameters Item are the same. For information on installing, viewing and editing assay parameters or for a detailed description of system procedures, refer to the Mindray Chemistry Analyzer Operator's Manual.

Calibration

1. For detailed calibration procedures, please refer to the Mindray Chemistry Analyzer Operator's Manual.

2. Calibration frequency: 8h by default.

Recalibration may be necessary when the followings occur:

- As ISE reagent lot changed;
- As the calibration factors are expired;
- As $\text{Na}^+/\text{K}^+/\text{Cl}^-$ /reference electrodes are replaced;
- As important components are replaced, such as sample tube, pump tube and tubing on the ISE module;
- As QC is out of control;
- After the ISE module and electrodes are maintained, including electrode cleaning and tube cleaning.

Each laboratory can set its own calibration procedures according to specific circumstances.

Quality control

Two levels of control material are recommended to analyze each

batch of samples. It is recommended to use the Mindray ClinChem Multi Control (Cat. No. 105-009119-00 and 105-009120-00) for internal serum/plasma quality control and MR Urine Quality Control (105-001691-00, Low: 1x100mL, High: 1x100mL) for internal urine quality control.

Controls provided by different manufacturers may be different. If a third-party QC is used, it is recommended to accumulate the target value and re-accumulate the target value when the lot of reagent and control is changed. The measured control value should be within the specified range. If the result is beyond the specified range, run the selected control at least once. If the result cannot be changed, it is recommended to replace the control or run it after recalibration. If necessary, take corresponding measures or contact the manufacturer.

Expected values

Sample Type		Range
Serum/Plasma	Adults	99-110 mmol/L ^{6,12,13,14}
Urine(24h)	Adults	110-250 mmol/24h ⁷

The expected value is provided from reference, and Mindray has verified it by at least 120 samples.

The urinary excretion of sodium, potassium and chloride varies significantly with dietary intake. The values given here are typical of people on an average diet.

Each laboratory should establish its own reference intervals based upon its local and population characteristics since expected values

may vary with geography, race, sex and age.

Performance Characteristics

■ Precision

Precision was determined by following CLSI Approved Guideline EP05-A3,⁸ control at two concentration levels and serum specimens at three concentration levels were assayed. Urine control at two levels are selected as urine specimen. Repeatability (n=20) was tested first, and then the specimen was assayed consecutively for 20 days. Each specimen was assayed 2 times per run and 2 runs per day (n=80). The precision data of controls and specimens on BS-600M are summarized below. Due to different specimen, instruments and operations, the data obtained in different laboratories may vary.

Serum:

Specimen Type (N=80)	Mean (mmol/L)	Repeatability		Within-Lab	
		SD mmol/L	CV %	SD mmol/L	CV %
Control1	76.9	0.3	0.3	0.8	1.1
Control2	105.7	0.4	0.3	0.3	0.3
Serum1	87.0	0.1	0.2	0.3	0.3
Serum2	107.5	0.2	0.1	0.3	0.3

Urine:

Specimen Type (N=80)	Mean (mmol/L)	Repeatability		Within-Lab	
		SD mmol/L	CV %	SD mmol/L	CV %
Control1	63.3	0.1	0.2	0.6	0.9
Control2	123.9	0.2	0.2	0.9	0.7

* Representative data or results in different instruments or laboratories may vary.

■ Measuring range

The Mindray chemistry analyzer provide the following linearity range:

Sample Type	Measuring range
Serum/Plasma	50-150
Urine	15-400

A high concentration Cl^- serum sample (approximately 160 mol/L) is mixed with a low concentration sample (<50 mmol/L) at different ratios, generating a series of dilutions. A high concentration Cl^- urine sample (approximately 440 mol/L) is mixed with a low concentration sample (<15 mmol/L) at different ratios, generating a series of dilutions. The Cl^- concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient $r \geq 0.990$.

■ Analytical Specificity

Interference studies were conducted using CLSI Approved Guideline EP07.⁹

The samples with different concentration interfering substances were prepared by addition of interferent to human sample pools, and if recovery is within $\pm 10\%$ of the corresponding control, it is considered as no significant interference.

No significant interference was observed when the following substances were tested for interference with this methodology. The

interference studies data of serum samples on BS-600M are summarized below.

Serum:

Interfering Substance	Interferent Concentration (mg/dL)	Analyte Concentration (mmol/L)	Relative Deviation (%) [*]
Hemoglobin	1000	95.59	-1.84
Unconjugated bilirubin	60	109.4	+2.19
Conjugated bilirubin	60	106.8	+5.11
Intralipid	2000	80.11	+4.04

^{}Representative data or results in different instruments or laboratories may vary.*

■ Method Comparison

Correlation studies were performed using CLSI Approved Guideline EP9-A3.¹⁰ The Mindray BS-600M system (y) was compared with Beckman AU5800 System using the same specimens. The statistical data obtained by linear regression are shown in the table below ^{*}:

The comparison data of Cl⁻ in serum specimens is as follows:

Regression Fit	Correlation Coefficient (r)	Sample (N)	Concentration Range (mg/L)
$y=1.022x-1.953$	0.962	120	77.3-125.1

The comparison data of Cl⁻ in urine specimens is as follows:

Regression Fit	Correlation	Sample (N)	Concentration
----------------	-------------	------------	---------------

Coefficient (r)		Range (mg/L)	
$y=0.971x+3.226$	0.990	101	15.9-313.4

* Representative data or results in different instruments or laboratories may vary.

Result interpretation

The results could be affected by drugs, disease, or endogenous substances. When the reaction curve is abnormal, it is recommended to retest and check the result.

Warnings and precautions

1. For in vitro diagnostic use only. For laboratory professional use. The installation and maintenance of electrodes must be operated by skilled/trained clinical professionals.
2. Please take the necessary precautions for handling all laboratory reagents.
3. Reliability of assay results cannot be guaranteed if the instructions in this package insert is not followed.
4. Please confirm the integrity of the package before use. Do not use the kits with damaged packages. The results can't be assured when stored at inappropriate condition.
5. Slope range of Cl^- electrode: 40-70 mV/decade. If the slope of the electrode exceeds the allowable range due to factors other than reagent, replace the electrode immediately.
6. Do not freeze the electrodes. Do not use electrodes and reagent that have expired. For the expiration date, see the package label. For the expiration date, see the package label. For electrode

replacement procedures, refer to the Mindray Chemistry Analyzer Operator's Manual. If the electrodes are stored refrigerated (2-8°C), they should be kept for a period of time and restored to room temperature before use.

7. After a large number of tests are finished, the proteins and lipids in the sample will adhere to the electrode surface, which may affect the measurement performance. It is recommended to clean the electrodes after finishing the ISE tests every day or maintain the electrodes before powering off the analyzer. For the specific cleaning steps, please refer to the Mindray Chemistry Analyzer Operator's Manual.
8. Biohazards will exist as electrodes may touch the patient's sample during the test, and the residual serum sample may contain virus and serum. Wear gloves when using the electrodes to avoid infection. If the sample touches the skin, wash the body area immediately with soap and water.
9. Before installing the electrodes, check visually whether the black sealing rubber ring of each electrode falls off, or it may cause poor sealing of the electrode tubes and affect the reliability of the results.
10. Before installing the electrodes, check visually if there are bubbles in the reference solution. Shake the electrodes gently to float up the bubbles when they exist.
11. Disposal of all samples and waste shall be in accordance with national and local laws and regulations.
12. Safety data sheet is available for professional user on request.

13. Install the electrodes according to the labels on the electrode compartment. Incorrect installation may result in calibration failure.
14. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration/controls do not meet the insert and/or the Mindray System criteria.
15. All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
16. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

References

1. Tietz NW, Pruden EL, Siggaard-Andersen O. Electrolytes. In: Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry. 2nd ed. Philadelphia: WB Saunders 1994: 1354-1374.
2. Eisenman G. Glass Electrodes for Hydrogen and Other Cations, Principle and Practice. New York: Marcel Dekker Inc.; 1967:2.
3. Kaplan L, Pesce A. Clinical Chemistry theory, analysis, and correlation St. Louis, Mosbby Co 1984:1061-1077.
4. Guder WG, Narayanan S, Wisser H, et al. List of analytes-preanalytical variables. Annex In: Samples: From the Patient to the Laboratory. Darmstadt, Germany: GIT Verlag; 1996: Annex 10-1, 20-1, 40-1.

5. US Pharmacopeial Convention, Inc. General notices. In: US Pharmacopeia National Formulary, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
6. WS/T404.3-2012 Reference Interval for Common Clinical Biochemistry Tests-Part 3: Serum Potassium, Sodium and Chloride.
7. Tietz NW. Fundamentals of Clinical Chemistry, 5th ed. Burtis CA, Ashwood ER, eds. WB Saunders Co 2001:970,1004,1009.
8. EP05-A3. Evaluation of Precision of Quantitative measurement procedures, 3rd.,2014.
9. EP07.Interference Testing in Clinical Chemistry, 3rd.,2018.
- 10.EP09-A3. Measurement procedure comparison and Bias Estimation Using Patient Samples, 3rd.,2013.
- 11.WHO/DIL/LAB/99.1 Rev.2,2002: Use of Anticoagulants in Diagnostic Laboratory Investigations& Stability of blood, plasma and serum samples. Geneva, 15 January 2002,pp41,48.
- 12.[U.S.A.] Carl ·A· Berty David ·E· Brens .Tietz Basis of Clinical Chemistry and Molecular Diagnostics [M]. Panbeshen. Chinese Medical Electronic Audio&Video Press, 2017.9: 409 -41.
- 13.Jia K, Zhang C, Huang X, Wang L, Hao X, Mu R, Pan B, Zhang J, Chen W, Xu N, Li G, Ma Y, Ma M, Guo W, Shang H. Reference Intervals of Serum Sodium, Potassium, and Chlorine in Chinese Han Population and Comparison of Two ISE Methods. J Clin Lab Anal. 2015 May;29(3):226-34. doi: 10.1002/jcla.21755. Epub 2014 May 5. PMID: 24799148; PMCID: PMC6807048.
- 14.Rustad P,Felding P,Franzson L et al. The Nordic Reference Interval Project 2000: recommended.

Graphical symbols

In vitro
diagnostic
medical device

Authorized
representative in
the European
Community

Consult
instructions for
use

European
Conformity

Batch code

Temperature limit

Manufacturer

Use-by date

Unique device
identifier

© 2022-2024 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech

Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998; **Fax:** +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH
(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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- Biochemistry Tests-Part 3: Serum Potassium, Sodium and Chloride.
7. Tietz NW. Fundamentals of Clinical Chemistry, 5th ed. Burtis CA, Ashwood ER, eds. WB Saunders Co 2001:970,1004,1009.
 8. EP05-A3. Evaluation of Precision of Quantitative measurement procedures, 3rd.,2014.
 9. EP07.Interference Testing in Clinical Chemistry, 3rd.,2018.
 - 10.EP09-A3. Measurement procedure comparison and Bias Estimation Using Patient Samples, 3rd.,2013.
 - 11.WHO/DIL/LAB/99.1 Rev.2,2002: Use of Anticoagulants in Diagnostic Laboratory Investigations& Stability of blood, plasma and serum samples. Geneva, 15 January 2002,pp41,48.
 - 12.[U.S.A.] Carl -A- Berty David -E- Brens .Tietz Basis of Clinical Chemistry and Molecular Diagnostics [M]. Panbeshen. Chinese Medical Electronic Audio&Video Press, 2017.9: 409 -41.
 - 13.Jia K, Zhang C, Huang X, Wang L, Hao X, Mu R, Pan B, Zhang J, Chen W, Xu N, Li G, Ma Y, Ma M, Guo W, Shang H. Reference Intervals of Serum Sodium, Potassium, and Chlorine in Chinese Han Population and Comparison of Two ISE Methods. J Clin Lab Anal. 2015 May;29(3):226-34. doi: 10.1002/jcla.21755. Epub 2014 May 5. PMID: 24799148; PMCID: PMC6807048.
 - 14.Rustad P,Felding P,Franzson L et al. The Nordic Reference Interval Project 2000: recommended.

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Referencias

1. Tietz NW, Pruden EL, Siggaard-Andersen O. Electrolytes. In: Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry. 2nd ed. Philadelphia: WB Saunders 1994: 1354-1374.
2. Eisenman G. Glass Electrodes for Hydrogen and Other Cations, Principle and Practice. New York: Marcel Dekker Inc.; 1967:2.
3. Kaplan L, Pesce A. Clinical Chemistry theory, analysis, and correlation St. Louis, Mosbby Co 1984:1061-1077.
4. Guder WG, Narayanan S, Wisser H, et al. List of analytes-preanalytical variables. Annex In: Samples: From the Patient to the Laboratory. Darmstadt, Germany: GIT Verlag; 1996: Annex 10-1, 20-1, 40-1.
5. US Pharmacopeial Convention, Inc. General notices. In: US Pharmacopeia National Formulary, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
6. WS/T404.3-2012 Reference Interval for Common Clinical Biochemistry Tests-Part 3: Serum Potassium, Sodium and Chloride.
7. Tietz NW. Fundamentals of Clinical Chemistry, 5th ed. Burtis CA, Ashwood ER, eds. WB Saunders Co 2001:970,1004,1009.
8. EP05-A3. Evaluation of Precision of Quantitative measurement procedures, 3rd.,2014.
9. EP07.Interference Testing in Clinical Chemistry, 3rd.,2018.
- 10.EP09-A3. Measurement procedure comparison and Bias

- Estimation Using Patient Samples, 3rd.,2013.
- 11.WHO/DIL/LAB/99.1 Rev.2,2002: Use of Anticoagulants in Diagnostic Laboratory Investigations& Stability of blood, plasma and serum samples. Geneva, 15 January 2002,pp41,48.
 - 12.[U.S.A.] Carl ·A· Berty David ·E· Brens ·Tietz Basis of Clinical Chemistry and Molecular Diagnostics [M]. Panbeshen. Chinese Medical Electronic Audio&Video Press, 2017.9: 409 -41.
 - 13.Jia K, Zhang C, Huang X, Wang L, Hao X, Mu R, Pan B, Zhang J, Chen W, Xu N, Li G, Ma Y, Ma M, Guo W, Shang H. Reference Intervals of Serum Sodium, Potassium, and Chlorine in Chinese Han Population and Comparison of Two ISE Methods. J Clin Lab Anal. 2015 May;29(3):226-34. doi: 10.1002/jcla.21755. Epub 2014 May 5. PMID: 24799148; PMCID: PMC6807048.
 - 14.Rustad P,Felding P,Franzson L et al. The Nordic Reference Interval Project 2000: recommended.

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Riferimenti

1. Tietz NW, Pruden EL, Siggaard-Andersen O. Electrolytes. In: Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry. 2nd ed.

-
- Philadelphia: WB Saunders 1994: 1354-1374.
2. Eisenman G. Glass Electrodes for Hydrogen and Other Cations, Principle and Practice. New York: Marcel Dekker Inc.; 1967:2.
 3. Kaplan L, Pesce A. Clinical Chemistry theory, analysis, and correlation St. Louis, Mosby Co 1984:1061-1077.
 4. Guder WG, Narayanan S, Wisser H, et al. List of analytes-preanalytical variables. Annex In: Samples: From the Patient to the Laboratory. Darmstadt, Germany: GIT Verlag; 1996: Annex 10-1, 20-1, 40-1.
 5. US Pharmacopeial Convention, Inc. General notices. In: US Pharmacopeia National Formulary, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
 6. WS/T404.3-2012 Reference Interval for Common Clinical Biochemistry Tests-Part 3: Serum Potassium, Sodium and Chloride.
 7. Tietz NW. Fundamentals of Clinical Chemistry, 5th ed. Burtis CA, Ashwood ER, eds. WB Saunders Co 2001:970,1004,1009.
 8. EP05-A3. Evaluation of Precision of Quantitative measurement procedures, 3rd.,2014.
 9. EP07.Interference Testing in Clinical Chemistry, 3rd.,2018.
 - 10.EP09-A3. Measurement procedure comparison and Bias Estimation Using Patient Samples, 3rd.,2013.
 - 11.WHO/DIL/LAB/99.1 Rev.2,2002: Use of Anticoagulants in Diagnostic Laboratory Investigations& Stability of blood, plasma and serum samples. Geneva, 15 January 2002,pp41,48.
 - 12.[U.S.A.] Carl ·A· Berty David ·E· Brens ·Tietz Basis of Clinical Chemistry and Molecular Diagnostics [M]. Panbeshen. Chinese
-

Medical Electronic Audio&Video Press, 2017.9: 409 -41.

- 13.Jia K, Zhang C, Huang X, Wang L, Hao X, Mu R, Pan B, Zhang J, Chen W, Xu N, Li G, Ma Y, Ma M, Guo W, Shang H. Reference Intervals of Serum Sodium, Potassium, and Chlorine in Chinese Han Population and Comparison of Two ISE Methods. J Clin Lab Anal. 2015 May;29(3):226-34. doi: 10.1002/jcla.21755. Epub 2014 May 5. PMID: 24799148; PMCID: PMC6807048.
- 14.Rustad P,Felding P,Franzson L et al. The Nordic Reference Interval Project 2000: recommended.

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Fabbricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Indirizzo: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Referanslar

- 1.Tietz NW, Pruden EL, Siggaard-Andersen O. Electrolytes. In: Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry. 2nd ed. Philadelphia: WB Saunders 1994: 1354-1374.
2. Eisenman G. Glass Electrodes for Hydrogen and Other Cations, Principle and Practice. New York: Marcel Dekker Inc.; 1967:2.
3. Kaplan L, Pesce A. Clinical Chemistry theory, analysis, and correlation St. Louis, Mosby Co 1984:1061-1077.
4. Guder WG, Narayanan S, Wisser H, et al. List of analytes-preanalytical variables. Annex In: Samples: From the Patient to the Laboratory. Darmstadt, Germany: GIT Verlag; 1996: Annex 10-1, 20-1, 40-1.
5. US Pharmacopeial Convention, Inc. General notices. In: US Pharmacopeia National Formulary, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
6. WS/T404.3-2012 Reference Interval for Common Clinical Biochemistry Tests-Part 3: Serum Potassium, Sodium and Chloride.
7. Tietz NW. Fundamentals of Clinical Chemistry, 5th ed. Burtis CA, Ashwood ER, eds. WB Saunders Co 2001:970,1004,1009.
8. EP05-A3. Evaluation of Precision of Quantitative measurement

- procedures, 3rd.,2014.
9. EP07.Interference Testing in Clinical Chemistry, 3rd.,2018.
 - 10.EP09-A3. Measurement procedure comparison and Bias Estimation Using Patient Samples, 3rd.,2013.
 - 11.WHO/DIL/LAB/99.1 Rev.2,2002: Use of Anticoagulants in Diagnostic Laboratory Investigations& Stability of blood, plasma and serum samples. Geneva, 15 January 2002,pp41,48.
 - 12.[U.S.A.] Carl ·A· Berty David ·E· Brens .Tietz Basis of Clinical Chemistry and Molecular Diagnostics [M]. Panbeshen. Chinese Medical Electronic Audio&Video Press, 2017.9: 409 -41.
 - 13.Jia K, Zhang C, Huang X, Wang L, Hao X, Mu R, Pan B, Zhang J, Chen W, Xu N, Li G, Ma Y, Ma M, Guo W, Shang H. Reference Intervals of Serum Sodium, Potassium, and Chlorine in Chinese Han Population and Comparison of Two ISE Methods. J Clin Lab Anal. 2015 May;29(3):226-34. doi: 10.1002/jcla.21755. Epub 2014 May 5. PMID: 24799148; PMCID: PMC6807048.
 - 14.Rustad P,Felding P,Franzson L et al. The Nordic Reference Interval Project 2000: recommended.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial
Park, Nanshan, Shenzhen, 518057, P. R. China

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Références

1. Tietz NW, Pruden EL, Siggaard-Andersen O. Electrolytes. In: Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry. 2nd ed. Philadelphia: WB Saunders 1994: 1354-1374.
2. Eisenman G. Glass Electrodes for Hydrogen and Other Cations, Principle and Practice. New York: Marcel Dekker Inc.; 1967:2.
3. Kaplan L, Pesce A. Clinical Chemistry theory, analysis, and correlation St. Louis, Mosby Co 1984:1061-1077.
4. Guder WG, Narayanan S, Wisser H, et al. List of analytes-preanalytical variables. Annex In: Samples: From the Patient to the Laboratory. Darmstadt, Germany: GIT Verlag; 1996: Annex 10-1, 20-1, 40-1.
5. US Pharmacopeial Convention, Inc. General notices. In: US Pharmacopeia National Formulary, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
6. WS/T404.3-2012 Reference Interval for Common Clinical Biochemistry Tests-Part 3: Serum Potassium, Sodium and Chloride.
7. Tietz NW. Fundamentals of Clinical Chemistry, 5th ed. Burtis CA, Ashwood ER, eds. WB Saunders Co 2001:970,1004,1009.
8. EP05-A3. Evaluation of Precision of Quantitative measurement procedures, 3rd., 2014.

9. EP07.Interference Testing in Clinical Chemistry, 3rd.,2018.
- 10.EP09-A3. Measurement procedure comparison and Bias Estimation Using Patient Samples, 3rd.,2013.
- 11.WHO/DIL/LAB/99.1 Rev.2,2002: Use of Anticoagulants in Diagnostic Laboratory Investigations& Stability of blood, plasma and serum samples. Geneva, 15 January 2002,pp41,48.
- 12.[U.S.A.] Carl -A- Berty David -E- Brens .Tietz Basis of Clinical Chemistry and Molecular Diagnostics [M]. Panbeshen. Chinese Medical Electronic Audio&Video Press, 2017.9: 409 -41.
- 13.Jia K, Zhang C, Huang X, Wang L, Hao X, Mu R, Pan B, Zhang J, Chen W, Xu N, Li G, Ma Y, Ma M, Guo W, Shang H. Reference Intervals of Serum Sodium, Potassium, and Chlorine in Chinese Han Population and Comparison of Two ISE Methods. J Clin Lab Anal. 2015 May;29(3):226-34. doi: 10.1002/jcla.21755. Epub 2014 May 5. PMID: 24799148; PMCID: PMC6807048.
- 14.Rustad P,Felding P,Franzson L et al. The Nordic Reference Interval Project 2000: recommended.

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Fabricant: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adresse: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Principle and Practice. New York: Marcel Dekker Inc.; 1967:2.

3. Kaplan L, Pesce A. Clinical Chemistry theory, analysis, and correlation St. Louis, Mosby Co 1984:1061-1077.
4. Guder WG, Narayanan S, Wisser H, et al. List of analytes-preanalytical variables. Annex In: Samples: From the Patient to the Laboratory. Darmstadt, Germany: GIT Verlag; 1996: Annex 10-1, 20-1, 40-1.
5. US Pharmacopeial Convention, Inc. General notices. In: US Pharmacopeia National Formulary, 1995 ed (USP 23/NF 18). Rockville, MD: The US Pharmacopeial Convention, Inc; 1994:11.
6. WS/T404.3-2012 Reference Interval for Common Clinical Biochemistry Tests-Part 3: Serum Potassium, Sodium and Chloride.
7. Tietz NW. Fundamentals of Clinical Chemistry, 5th ed. Burtis CA, Ashwood ER, eds. WB Saunders Co 2001:970,1004,1009.
8. EP05-A3. Evaluation of Precision of Quantitative measurement procedures, 3rd.,2014.
9. EP07.Interference Testing in Clinical Chemistry, 3rd.,2018.
- 10.EP09-A3. Measurement procedure comparison and Bias Estimation Using Patient Samples, 3rd.,2013.
- 11.WHO/DIL/LAB/99.1 Rev.2,2002: Use of Anticoagulants in Diagnostic Laboratory Investigations& Stability of blood, plasma and serum samples. Geneva, 15 January 2002,pp41,48.
- 12.[U.S.A.] Carl -A- Berty David -E- Brens .Tietz Basis of Clinical Chemistry and Molecular Diagnostics [M]. Panbeshen. Chinese Medical Electronic Audio&Video Press, 2017.9: 409 -41.
- 13.Jia K, Zhang C, Huang X, Wang L, Hao X, Mu R, Pan B, Zhang J,

Chen W, Xu N, Li G, Ma Y, Ma M, Guo W, Shang H. Reference Intervals of Serum Sodium, Potassium, and Chlorine in Chinese Han Population and Comparison of Two ISE Methods. J Clin Lab Anal. 2015 May;29(3):226-34. doi: 10.1002/jcla.21755. Epub 2014 May 5. PMID: 24799148; PMCID: PMC6807048.

14.Rustad P,Felding P,Franzson L et al. The Nordic Reference Interval Project 2000: recommended.

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Producent: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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