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HEMOGLOBIN A1C CONTROL (NORMAL, ELEVATED)

NORMAL	ELEVATED
COD 18001 1 x 0.5 mL	COD 18002 1 x 0.5 mL
STORE AT 2-8°C	
Only for <i>in vitro</i> use in the clinical laboratory	

INTENDED USE

The Hemoglobin A_{1c} Control intended for quality control purposes only and is supplied with suggested intervals of acceptable values.

COMPOSITION

Hemoglobin A_{1c} Control. For 0.5 mL. Lyophilized hemolysates from human origin, suitable for the quality control of the clinical laboratories, and without preservatives which might interfere with the tests.

Components from human origin have been tested and found to be negative for the presence of antibodies anti-HIV and anti-HCV, as well as for HBs antigen. However, they should be handled cautiously as potentially infectious.

PREPARATION AND USE

1. Open a vial very carefully, avoiding any loss of the lyophilized material.
2. Pipette exactly 0.5 mL of distilled water into the vial. The component values depend on the accuracy of this reconstitution step.
3. Close the vial with the stopper and let stand for 20 minutes at room temperature.
4. Swirl gently, avoiding the formation of foam, to ensure complete dissolution of contents.
5. Analyze the reconstituted material in the same way as the patient samples.

STORAGE AND STABILITY

Store at 2-8°C.

The lyophilized Control is stable until the expiration date given in the label.

The components of the reconstituted material are stable for at least 7 days at 2-8°C. For longer periods, keep frozen. Avoid repeated freeze-thaw cycles.

SPECIFIED VALUES

Hemoglobin A_{1c} concentration is given on the vial label, expressed in mmol/L and %. The value is traceable to the reference method as described by IFCC. Traceability of the results can be assured only if the BioSystems reagents and recommended measurement procedures are used.

The suggested intervals of acceptable values have been calculated from previous experience in interlaboratory variability and are given for orientation only. Each laboratory should establish its own precision parameters.

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08/2024

M180011-08

The lateral lines mark the modifications in the current version.

08/2024

BioSystems S.A. Costa Brava, 30. 08030 Barcelona (Spain)
www.biosystems.global

HEMOGLOBIN A1C CONTROL (NORMAL, ELEVATED)

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