

Detro® Cid Enzym

Enzymatic Cleaner and Disinfectant for Medical Instruments and Endoscopes



USAGE AREAS:

- Used for cleaning and disinfection of surgical instruments, felxibel and rigid endoscopes. It is compatible with medical devices/equipments made of polystyrene, stainless steel, PVC, acrylic material, polycarbon, ceramic, ABS, silicone, HDPE, rubber and glass.

PRODUCT PROPERTIES:

- Concentrated product.
- Disinfection and cleaning in one step.
- Effective cleaner thanks to enzymes, cleans organic dirt like blood, protein and oil.
- Suitable for ultrasonic and manual usage.
- Aldehyde and phenol free.
- Prevents corrosion of tools thanks to corrosion inhibitors, compatible with endoscope brands.

COMPOSITION:

- Enzymes (amylase, lipase, protease), didecyldimethylammonium chloride, benzalconium chloride, dodecylamine, corrosion inhibitors, auxiliary substances.

MICROBIOLOGICAL ACTIVITY:

- Bactericidal, fungicidal, virucida, tuberculocidal in 15 minutes at %0,5 concentration.

PACKAGE:

QUANTITY	PACKING TYPE	BARCODE
2L	6 x BOTTLE	8680097100719
5L	4 x DRUM	8680097100511



USAGE:

- Detro® Cid Enzym is diluted at proper ratio and filled into a clean and closed container.
- The instruments are immersed thoroughly so that all the surfaces and cavities are covered by the solution.
- The container is closed.
- The contact time is completed.
- Then the instruments are rinsed with sterile or distilled water.
- The ready to use solution in the container can be used 15 days, the activity should be controlled by test strips, the solution should be renewed when the strips doesn't give the active colour exisiting on the strip tubes.



USAGE AREAS:

- Used for cleaning and disinfection of surgical instruments, felxibel and rigid endoscopes. It is compatible with medical devices/equipments made of polystyrene, stainless steel, PVC, acrylic material, polycarbon, ceramic, ABS, silicone, HDPE, rubber and glass.

PRODUCT PROPERTIES:

- Concentrated product.
- Disinfection and cleaning in one step.
- Suitable for ultrasonic and manual usage.
- Aldehyde and phenol free.
- Prevents corrosion of tools thanks to corrosion inhibitors, compatible with endoscope brands.

COMPOSITION:

- Bis (3-aminopropyl) dodecylamine, didecyldimethylammonium chloride, benzalconium chloride, corrosion inhibitors, auxiliary substances

MICROBIOLOGICAL ACTIVITY:

- Bactericidal, fungicidal, virucida, tuberculocidal in 15 minutes at %0,5 concentration.

PACKAGE:

QUANTITY	PACKING TYPE	BARCODE
1L	12 x BOTTLE	8680097103192
2L	6 x BOTTLE	8680097101839
5L	4 x DRUM	8680097100702



USAGE:

- Diluted at proper ratio and filled into a clean and closed container.
- The instruments are immersed thoroughly so that all the surfaces and cavities are covered by the solution.
- The container is closed.
- The contact time is completed.
- Then the instruments are rinsed with sterile or distilled water.
- The ready to use solution in the container can be used 15 days, the activity should be controlled by test strips, the solution should be renewed when the strips doesn't give the active colour exisiting on the strip tubes.