EIAgen

HCV Ab (v.4) Kit

REF 071067

 96

REF 071064

 192

REF 071068

 480

IVD

CE 0459

This package insert must be read carefully before product use.

Package insert instructions must be carefully followed.

Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package.

Manufacturer:

Adaltis S.r.l

Via Durini, 27

20122 Milano (Italy)

Tel. +39-0774-5791 - Fax +39-0774-353085

www.adaltis.net

en

SYMBOLS USED ON LABELS

English EN							
	In Vitro Diagnostic Medical Device	Catalogue Number	Lot Number	Attention, See Instructions For Use	Temperature Limitation	Use By	Number of Test
	Manufacturer	Keep away from Sunlight	Date of Manufacture	Microplate	Positive Control	Negative Control	Calibrator
	Conjugate	Sample Diluent	Substrate TMB	Stop Solution (H ₂ SO ₄ 0,3 M)	Wash Buffer Conc. (20x)	Assay Diluent	Reconstitute with x mL
	Biological Risk	Danger	Warning				

Attention:

Stop Solution classified as: Skin Corr. 1A

- **Signal word:**
Danger
- **Hazard-determining components of labelling:**
Sulphuric acid
- **Hazard statements:**
H314 Causes severe skin burns and eye damage.
- **Precautionary statements:**
P260 Do not breathe dust/fume/gas/mist/vapours/spray.
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER/doctor.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Attention:

Negative Control, Positive Control, Calibrator, Conjugate, Sample Diluent, Assay Diluent and Wash Buffer Concentrate 20X classified as: Skin Sens. 1

- **Signal word:**
Warning
- **Hazard-determining components of labelling:**
Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)
- **Hazard statements:**
H317 May cause an allergic skin reaction.
- **Precautionary statements:**
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P321 Specific treatment (see on this label).
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P302+P352 IF ON SKIN: Wash with plenty of water.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Refer to www.adaltis.net for the Safety Data Sheets.

A. INTENDED USE

Fourth generation Enzyme ImmunoAssay (ELISA) for the determination of antibodies to Hepatitis C Virus in human plasma (EDTA, Heparin and Citrate) and sera. The kit may be used for the screening of blood units of HCV-infected patients.

For "in vitro" diagnostic use only.

B. INTRODUCTION

The World Health Organization (WHO) define Hepatitis C infection as follows:

"Hepatitis C is a viral infection of the liver which had been referred to as parenterally transmitted "non A, non B hepatitis" until identification of the causative agent in 1989. The discovery and characterization of the hepatitis C virus (HCV) led to the understanding of its primary role in post transfusion hepatitis and its tendency to induce persistent infection.

HCV is a major cause of acute hepatitis and chronic liver disease, including cirrhosis and liver cancer. Globally, an estimated 170 million persons are chronically infected with HCV and 3 to 4 million persons are newly infected each year. HCV is spread primarily by direct contact with human blood. The major causes of HCV infection worldwide are use of unscreened blood transfusions, and re-use of needles and syringes that have not been adequately sterilized. No vaccine is currently available to prevent hepatitis C and treatment for chronic hepatitis C is too costly for most persons in developing countries to afford. Thus, from a global perspective, the greatest impact on hepatitis C disease burden will likely be achieved by focusing efforts on reducing the risk of HCV transmission from nosocomial exposures (e.g. blood transfusions, unsafe injection practices) and high-risk behaviours (e.g. injection drug use).

Hepatitis C virus (HCV) is one of the viruses (A, B, C, D, and E), which together account for the vast majority of cases of viral hepatitis. It is an enveloped RNA virus in the flaviviridae family which appears to have a narrow host range. Humans and chimpanzees are the only known species susceptible to infection, with both species developing similar disease. An important feature of the virus is the relative mutability of its genome, which in turn is probably related to the high propensity (80%) of inducing chronic infection. HCV is clustered into several distinct genotypes which may be important in determining the severity of the disease and the response to treatment.

The incubation period of HCV infection before the onset of clinical symptoms ranges from 15 to 150 days. In acute infections, the most common symptoms are fatigue and jaundice; however, the majority of cases (between 60% and 70%), even those that develop chronic infection, are asymptomatic. About 80% of newly infected patients progress to develop chronic infection. Cirrhosis develops in about 10% to 20% of persons with chronic infection, and liver cancer develops in 1% to 5% of persons with chronic infection over a period of 20 to 30 years. Most patients suffering from liver cancer who do

not have hepatitis B virus infection have evidence of HCV infection. The mechanisms by which HCV infection leads to liver cancer are still unclear. Hepatitis C also exacerbates the severity of underlying liver disease when it coexists with other hepatic conditions. In particular, liver disease progresses more rapidly among persons with alcoholic liver disease and HCV infection. HCV is spread primarily by direct contact with human blood. Transmission through blood transfusions that are not screened for HCV infection, through the reuse of inadequately sterilized needles, syringes or other medical equipment, or through needle-sharing among drug-users, is well documented. Sexual and perinatal transmission may also occur, although less frequently. Other modes of transmission such as social, cultural, and behavioural practices using percutaneous procedures (e.g. ear and body piercing, circumcision, tattooing) can occur if inadequately sterilized equipment is used. HCV is not spread by sneezing, hugging, coughing, food or water, sharing eating utensils, or casual contact.

In both developed and developing countries, high risk groups include injecting drug users, recipients of unscreened blood, haemophiliacs, dialysis patients and persons with multiple sex partners who engage in unprotected sex. In developed countries, it is estimated that 90% of persons with chronic HCV infection are current and former injecting drug users and those with a history of transfusion of unscreened blood or blood products. In many developing countries, where unscreened blood and blood products are still being used, the major means of transmission are unsterilized injection equipment and unscreened blood transfusions. In addition, people who use traditional scarification and circumcision practices are at risk if they use or re-use unsterilized tools.

WHO estimates that about 170 million people, 3% of the world's population, are infected with HCV and are at risk of developing liver cirrhosis and/or liver cancer. The prevalence of HCV infection in some countries in Africa, the Eastern Mediterranean, South-East Asia and the Western Pacific (when prevalence data are available) is high compared to some countries in North America and Europe.

Diagnostic tests for HCV are used to prevent infection through screening of donor blood and plasma, to establish the clinical diagnosis and to make better decisions regarding medical management of a patient. Diagnostic tests commercially available today are based on Enzyme immunoassay (EIA) for the detection of HCV specific antibodies. EIAs can detect more than 95% of chronically infected patients but can detect only 50% to 70% of acute infections. A recombinant immunoblot assay (RIBA) that identifies antibodies which react with individual HCV antigens is often used as a supplemental test for confirmation of a positive EIA result. Testing for HCV circulating by amplification tests RNA (e.g. polymerase chain reaction or PCR, branched DNA assay) is also being utilized for confirmation of serological results as well as for assessing the effectiveness of antiviral therapy. A positive result indicates the presence of active infection and a potential for spread of the infection and or/the development of chronic liver disease.

Antiviral drugs such as interferon taken alone or in combination with ribavirin, can be used for the treatment of persons with chronic hepatitis C, but the cost of treatment is very high. Treatment with interferon alone is effective in about 10% to 20% of patients. Interferon combined with ribavirin is effective in about 30% to 50% of patients. Ribavirin does not appear to be effective when used alone.

There is no vaccine against HCV. Research is in progress but the high mutability of the HCV genome complicates vaccine development. Lack of knowledge of any protective immune response following HCV infection also impedes vaccine research. It is not known whether the immune system is able to eliminate the virus.

Some studies, however, have shown the presence of virus neutralizing antibodies in patients with HCV infection. In the absence of a vaccine, all precautions to prevent infection must be taken including (a) screening and testing of blood and organ donors; (b) Virus inactivation of plasma derived products; (c) implementation and maintenance of infection control practices in health care settings, including appropriate sterilization of medical and dental equipment; (d) promotion of behaviour change among the general public and health care workers to reduce overuse of injections and to use safe injection practices; and (e) Risk reduction counselling for persons with high-risk drug and sexual practices. “

The genome encodes for structural components, a nucleocapsid protein and two envelope glycoproteins, and functional constituents involved in the virus replication and protein processing. The nucleocapsid-encoding region seems to be the most conservative among the isolates obtained all over the world.

C. PRINCIPLE OF THE TEST

Microplates are coated with HCV-specific antigens derived from “core” and “ns” regions encoding for conservative and immunodominant antigenic determinants (Core peptide, recombinant NS3, NS4 and NS5 peptides).

The solid phase is first treated with the diluted sample and HCV Ab are captured, if present, by the antigens. After washing out all the other components of the sample, in the 2nd incubation bound HCV antibodies, IgG and IgM as well, are detected by the addition of polyclonal specific anti hlgG&M antibodies, labelled with peroxidase (HRP).

The enzyme captured on the solid phase, acting on the substrate TMB mixture, generates an optical signal that is proportional to the amount of anti HCV antibodies present in the sample. A cut-off value let optical densities be interpreted into HCV antibody negative and positive results.

D. COMPONENTS

The kit contains reagents for 96 tests (code 071067), 192 tests (code 071064), or 480 tests (code 071068).

Microplate	1
Negative Control	1x4 mL/vial
Positive Control	1x2 mL/vial
Calibrator	2 vials
Wash Buffer Concentrate 20x	1x50 mL/vial
Conjugate	1x16 mL/vial
Sample Diluent	1x50 mL/vial
Substrate TMB	1x16 mL/vial
Stop Solution	1x15 mL/vial
Assay Diluent	1x8 mL/vial
Plate sealing foils	2
Number of tests	96
Code	071067

Microplate	2
Negative Control	2x4 mL/vial
Positive Control	1x4 mL/vial
Calibrator	3 vials
Wash Buffer Concentrate 20x	2x50 mL/vials
Conjugate	2x16 mL/vials
Sample Diluent	2x50 mL/vials
Substrate TMB	2x16 mL/vials
Stop Solution	2x15 mL/vial
Assay Diluent	2x8 mL/vial
Plate sealing foils	4
Number of tests	192
Code	071064

Microplate	5
Negative Control	1x20 mL/vial
Positive Control	1x10 mL/vial
Calibrator	7 vials
Wash Buffer Concentrate 20x	5x50 mL/vials
Conjugate	2x40 mL/vial
Sample Diluent	5x50 mL/vials
Substrate TMB	2x40 mL/vials
Stop Solution	2x40 mL/vial
Assay Diluent	1x40 mL/vial
Plate sealing foils	10
Number of tests	480
Code	071068

1. Microplate

12 strips of 8 microwells coated with Core peptide, recombinant NS3, NS4 and NS5 peptides. Plates are sealed into a aluminium pouch with desiccant. Bring the microplate to room temperature (18...24°C) before opening the bag. Unused strips have to be returned into the pouch and the pouch has to be sealed and stored back to 2...8°C, in presence of the desiccant.

2. Negative Control

Ready to use control. It contains 10 mM Na-citrate buffer pH 6.0 ± 0.1, 2% casein as protein base and 0.1% Proclin 150 as preservative. The negative control is olive green colour coded.

3. Positive Control

Ready to use control. It contains 1% goat serum proteins, human antibodies positive to HCV, 10 mM Na-citrate buffer pH 6.0 ± 0.1, 0.5% Tween 20, 0.09% Na-azide and 0.1% Proclin 150 as preservatives. The Positive Control is dark green colour coded.

Important Note: The absence of viable pathogens in the Positive Control can not be fully ensured, and therefore,

the reagent should be handled as potentially biohazardous, in accordance with good laboratory practices.

4. Calibrator

Lyophilized calibrator. To be dissolved with the volume of EIA grade water reported on the label. It contains foetal bovine serum proteins, human antibodies to HCV whose content is calibrated on the NIBSC Working Standard code 06/188-006, 10 mM Na-citrate buffer pH 6.0 ± 0.1 , 0.3 mg/mL gentamicin sulfate and 0.1% Proclin 150 as preservatives.

Important Note: *The absence of viable pathogens in the Calibrator can not be fully ensured, and therefore, the reagent should be handled as potentially biohazardous, in accordance with good laboratory practices.*

Note: *The volume necessary to dissolve the content of the vial may vary from lot to lot. Please use the right volume reported on the label.*

5. Wash Buffer Concentrate 20x

20x concentrated solution. Once diluted, the wash solution (wash buffer diluted) contains 10 mM phosphate buffer pH 7.0 ± 0.2 , 0.05% Tween 20 and 0.05% Proclin 150.

Once diluted, the wash solution is stable for 1 week at 2...8°C.

6. Conjugate

Ready to use and red colour coded reagent. It contains Horseradish Peroxidase conjugated goat polyclonal antibodies to human IgG and IgM, 5% BSA, 10 mM Citrate buffer pH 6.4 ± 0.1 , 0.1% Proclin and 0.05% Tween 20 as preservatives.

7. Substrate TMB

Ready-to-use component. It contains 50 mM citrate-phosphate buffer pH 3.5-3.8, 4% dimethylsulphoxide, 0.03% tetra-methyl-benzidine or TMB and 0.02% hydrogen peroxide or H₂O₂. Mix gently before use.

Note: *To be stored protected from light as sensitive to strong illumination.*

8. Assay diluent

Ready-to-use component. It contains goat serum, 10 mM tris buffered solution pH 8.0 ± 0.1 containing 0.1% Proclin 150 and 0.09% Na azide for the pre-treatment of samples and controls in the plate, blocking interference.

9. Stop Solution

Ready-to-use component.

It contains 0.3 M H₂SO₄ solution. Mix gently before use.

10. Sample Diluent

Ready-to-use component and dark green colour coded. It contains 1% casein, 10 mM Na-citrate buffer pH 6.0 ± 0.1 and 0.1% Proclin 150 as preservative.

To be used to dilute the sample.

Note: *The diluent changes colour from olive green to dark bluish green in the presence of sample.*

E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes (200 μ L and 10 μ L) and disposable plastic tips.
2. EIA grade water (bidistilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator capable to provide a temperature of +37°C.
6. Calibrated ELISA microwell reader with 450nm (reading) and possibly with 620-630nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory. This package insert must be read carefully before product use.
2. Read carefully the Safety Data Sheet (SDS) before product use.
3. When the kit is used for the screening of blood units and blood components, it has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
4. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
5. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
6. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Substrate (TMB) from strong light and avoid vibration of the bench surface where the test is undertaken.
7. Upon receipt, store the kit at 2...8°C into a temperature controlled refrigerator or cold room.
8. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
9. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
10. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
11. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.

12. Do not use the kit after the expiration date stated on the external container and internal (vials) labels.
13. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
14. The use of disposable plastic-ware is recommended in the preparation of the liquid components or in transferring components into automated workstations, in order to avoid cross contamination.
15. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated before waste. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min.
16. Accidental spills from samples and operations have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
17. The Stop Solution contains 0.3M sulphuric acid. Avoid contact with skin and eyes. In the event of contact, rinse immediately with plenty of water.
18. Dispose of reagent solutions containing sodium azide as preservatives according to all local, state and national regulations. To dispose of reagents containing azide, flush away using copious amounts of water. Dispose with caution as sodium azide can form explosive compounds on prolonged contact with lead or copper piping.
19. Do not smoke, eat, drink or apply cosmetics in areas in which specimens or kit reagents are handled
20. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.
21. Do not pipette by mouth.

G. SPECIMEN: PREPARATION AND RECCOMANDATIONS

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.
2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labeling and electronic reading is strongly recommended.

4. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
5. Sera and plasma can be stored at +2...8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.

H. PREPARATION OF COMPONENTS AND WARNINGS

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 1 re-use of the device and up to 6 months.

1. Microplates:

Allow the microplate to reach room temperature (about 1 hr) before opening the pouch. Check that the desiccant has not turned to dark green, indicating a defect in manufacturing. In this case, call Adaltis's customer service.

Unused strips have to be placed back inside the aluminium pouch, with the desiccant supplied, firmly zipped and stored at +2...8°C. After first opening, remaining strips are stable until the humidity indicator inside the desiccant bag turns from yellow to green.

2. Negative Control:

Ready to use. Mix well on vortex before use.

3. Positive Control:

Ready to use. Mix well on vortex before use. Handle this component as potentially infectious.

4. Calibrator:

Dissolve carefully the content of the lyophilised vial with the volume of EIA grade water reported on its label. Mix well on vortex before use.

Handle this component as potentially infectious.

Note: When dissolved the Calibrator is not stable. Store in aliquots at -20°C.

5. Wash Buffer Concentrate 20x (vial of 50 mL):

The whole content of the 20x concentrated solution has to be diluted with EIA grade water up to 1000 mL (the volume is reported on the label) and mixed gently end-over-end before use. As some salt crystals may be present into the vial, take care to dissolve all the content when preparing the solution. During preparation avoid foaming as the presence of bubbles could impact on the efficiency of the washing cycles.

Note: Once diluted, the wash solution is stable for 1 week at +2...8°C.

6. Conjugate:

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes. If this component

has to be transferred use only plastic, possibly sterile disposable containers.

7. Substrate TMB:

Ready to use. Mix well on vortex before use.

Avoid contamination of the liquid with oxidizing chemicals, air-driven dust or microbes. Do not expose to strong illumination, oxidizing agents and metallic surfaces. If this component has to be transferred use only plastic, and if possible, sterile disposable container.

8. Assay Diluent:

Ready to use. Mix well on vortex before use.

9. Stop Solution:

Ready to use. Mix well on vortex before use.

10. Sample Diluent:

Ready to use. Mix well on vortex before use.

I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample or the components of the kit. They should also be regularly maintained in order to show a precision of 1% and a trueness of $\pm 2\%$. Decontamination of spills or residues of kit components should also be carried out regularly.
2. The ELISA incubator has to be set at $+37^{\circ}\text{C}$ (tolerance of $\pm 0.5^{\circ}\text{C}$) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must be carefully validated and correctly optimized. Usually 4-5 washing cycles (aspiration + dispensation of $350\ \mu\text{L}/\text{well}$ of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected. A soaking time of 20-30 seconds between cycles is suggested. In order to set correctly their number, it is recommended to run an assay with the kit controls/calibrator and well-characterized negative and positive reference samples, and check to match the values reported below in sections O "Internal Quality Control". Regular calibration of the volumes delivered and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.
4. Incubation times have a tolerance of $\pm 5\%$.
 - ✓ Short Incubation Method (for 1st / 2nd incubation tolerance between 43 min to 47 min; for 3rd incubation tolerance between 14 and 16 min).
 - ✓ Standard Incubation Method (for 1st incubation tolerance between 57 min to 63 min; for 2nd and 3rd incubation tolerance between 29 and 31 min).
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and ideally with a second

filter (620-630nm) for blanking purposes. Its standard performances should be (a) bandwidth $\leq 10\ \text{nm}$; (b) absorbance range from 0 to ≥ 2.0 ; (c) linearity to ≥ 2.0 ; repeatability $\geq 1\%$. Blanking is carried out on the well identified in section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer's instructions.

6. When using an ELISA automated work station, all critical steps (dispensation, incubation, washing, reading, shaking, data handling) have to be carefully set, calibrated, controlled and regularly serviced in order to match the values reported in sections O "Internal Quality Control". The assay protocol has to be installed in the operating system of the unit and validated as for the washer and the reader. In addition, the liquid handling part of the station (dispensation and washing) has to be validated and correctly set. Particular attention must be paid to avoid carry over by the needles used for dispensing samples and for washing. This must be studied and controlled to minimize the possibility of contamination of adjacent wells due to strongly reactive samples, leading to false positive results. The use of ELISA automated work stations is recommended for blood screening and when the number of samples to be tested exceed 20-30 units per run.

L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label of the kit box. Do not use if expired.
2. Check that the liquid components are not contaminated by naked-eye visible particles or aggregates. Check that the Substrate TMB is colorless or pale blue by aspirating a small volume of it with a sterile transparent plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminum pouch, containing the microplate, is not punctured or damaged.
3. Dilute all the content of the 20x concentrated Wash Buffer as described above.
4. Dissolve the Calibrator as described above.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix as described.
6. Set the ELISA incubator at $+37^{\circ}\text{C}$ and prepare the ELISA washer by priming with the diluted wash solution, according to the manufacturers instructions. Set the right number of washing cycles as found in section I.3.
7. Check that the ELISA reader has been turned on at least 20 minutes before reading.
8. If using an automated workstation, turn it on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipment is available and ready to use.
11. In case of problems, do not proceed further with the test and advise the supervisor.

M. ASSAY PROCEDURE

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

The Assay Procedure can be done with two times incubation procedures. Choose the one that is required by your regulation:

1. Standard Incubation (1st incubation 60 minutes, 2nd and 3rd incubation 30 minutes)
2. Short Incubation (1st and 2nd incubation 45 minutes, 3rd incubation 15 minutes)

1. Standard Incubation - Manual assay:

1. Place the required number of microwells in the microwell holder. Leave the 1st well empty for the operation of blanking.
2. Dispense 200 µL of Negative Control in triplicate, 200 µL Calibrator in duplicate and 200 µL Positive Control in single in proper wells. Do not dilute Controls and Calibrator as they are pre-diluted, ready to use!
3. Add 200 µL of Sample Diluent to all the sample wells; then dispense 10 µL sample in each properly identified well. Mix gently the plate, avoiding overflowing and contaminating adjacent wells, in order to fully disperse the sample into its diluent.

Important note: Check that the colour of the Sample Diluent, upon addition of the sample, changes from light green to dark bluish green, monitoring that the sample has been really added.

4. Dispense 50 µL Assay Diluent into all the controls/calibrator and sample wells. Check that the color of samples has turned to dark blue.
5. Incubate the microplate for **60 min at +37°C**.
Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.
6. Wash the microplate with an automatic washer by delivering and aspirating 350 µL/well of diluted wash solution as reported in section I.3.
7. Pipette 100 µL Enzyme Conjugate into each well, except the 1st blanking well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

8. Incubate the microplate for **30 min at +37°C**.
9. Wash microwells as reported in section I.3.
10. Pipette 100 µL Substrate TMB mixture into each well, the blank well included. Then incubate the microplate at room temperature (18...24°C) for **30 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

11. Pipette 100 µL Stop Solution into all the wells using the same pipetting sequence as in step 10 to stop the enzymatic reaction. Addition of stop solution will

turn the positive control and positive samples from blue to yellow.

12. Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and possibly at 620-630nm (background subtraction), blanking the instrument on A1.

Important notes:

1. If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
2. Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the substrate can occur leading to a higher background.
3. Shaking at 350 ±150 rpm during incubation has been proved to increase the sensitivity of the assay of about 20%.

2. Short Incubation - Manual assay:

1. Place the required number of Microwells in the microwell holder. Leave the 1st well empty for the operation of blanking.
2. Dispense 200 µL of Negative Control in triplicate, 200 µL Calibrator in duplicate and 200 µL Positive Control in single in proper wells. Do not dilute Controls and Calibrator as they are pre-diluted, ready to use!
3. Add 200 µL of Sample Diluent to all the sample wells; then dispense 10 µL sample in each properly identified well. Mix gently the plate, avoiding overflowing and contaminating adjacent wells, in order to fully disperse the sample into its diluent.

Important note: Check that the colour of the Sample Diluent, upon addition of the sample, changes from light green to dark bluish green, monitoring that the sample has been really added.

4. Dispense 50 µL Assay Diluent into all the controls/calibrator and sample wells. Check that the color of samples has turned to dark blue.
5. Incubate the microplate for **45 min at +37°C**.
Important note: Strips have to be sealed with the adhesive sealing foil, supplied, only when the test is carried out manually. Do not cover strips when using ELISA automatic instruments.
6. Wash the microplate with an automatic washer by delivering and aspirating 350 µL/well of diluted washing solution as reported in section I.3.
7. Pipette 100 µL Enzyme Conjugate into each well, except the 1st blanking well, and cover with the sealer. Check that this red coloured component has been dispensed in all the wells, except A1.

Important note: Be careful not to touch the plastic inner surface of the well with the tip filled with the Enzyme Conjugate. Contamination might occur.

8. Incubate the microplate for **45 min at +37°C**.
9. Wash microwells as reported in section I.3.
10. Pipette 100 µL Substrate TMB mixture into each well, the blank well included. Then incubate the

microplate at room temperature (18...24°C) for **15 minutes**.

Important note: Do not expose to strong direct illumination. High background might be generated.

11. Pipette 100 µL Stop Solution into all the wells using the same pipetting sequence as in step 10 to stop the enzymatic reaction. Addition of stop solution will turn the positive control and positive samples from blue to yellow.
12. Measure the colour intensity of the solution in each well, as described in section I.5, at 450nm filter (reading) and possibly at 620-630nm (background subtraction), blanking the instrument on A1.

Important notes:

1. If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
2. Reading has should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the substrate can occur leading to a higher background.
3. Shaking at 350 ±150 rpm during incubation has been proved to increase the sensitivity of the assay of about 20%.

N. ASSAY SCHEME

Method	Operations (Standard Incubation)	Operations (Short Incubation)
Controls & Calibrator	200 µL	200 µL
Sample Diluent & Sample	200 µL diluent.+ 10 µL sample	200 µL diluent.+ 10 µL sample
Assay Diluent	50 µL	50 µL
1st incubation	60 min (± 3)	45 min (± 2)
Temperature	+37°C	+37°C
Wash step	4-5 cycles	4-5 cycles
Enzyme Conjugate	100 µL	100 µL
2nd incubation	30 min (± 1)	45 min (± 2)
Temperature	+37°C	+37°C
Wash step	4-5 cycles	4-5 cycles
Substrate TMB	100 µL	100 µL
3rd incubation	30 min (± 1)	15 min (± 1)
Temperature	Room temperature (18...24°C)	Room temperature (18...24°C)
Stop Solution	100 µL	100 µL
Reading OD	450/620nm	450/620nm

An example of dispensation scheme is reported below (valid for both incubation time procedures):

Microplate

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	S2										
B	NC	S3										
C	NC	S4										
D	NC	S5										
E	CAL	S6										
F	CAL	S7										
G	PC	S8										
H	S1	S9										

Legenda: BLK = Blank NC = Negative Control
CAL = Calibrator PC = Positive Control S = Sample

O. INTERNAL QUALITY CONTROL

A check is carried out on the controls and the calibrator any time the kit is used in order to verify whether their OD450/620nm values are as expected and reported in the table below.

Check	Requirements
Blank well	< 0.100 OD 450/620nm value
Negative Control (NC)	< 0.050 mean OD450/620nm after blanking
Calibrator	S/Co >1.1
Positive Control	>1.000 OD450/620nm value

If the results of the test match the requirements stated above, proceed to the next section. If they do not, do not proceed any further and perform the following checks:

Problems	Check
Blank well > 0.100 OD450nm	1. that the Substrate solution has not become contaminated during the assay
Negative Control (NC) > 0.050 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of the negative one); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to spills of positive samples or of the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Calibrator S/Co<1.1	1. that the procedure has been correctly performed; 2. that no mistake has occurred during its distribution (ex.: dispensation of negative control instead of calibrator) 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the calibrator has occurred.
Positive Control <1.000 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control. In this case, the negative control will have an OD450nm value > 0.150,too. 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

P. RESULTS

The tests results are calculated by means of a cut-off value determined with the following formula:

$$\text{Cut-Off (Co)} = \text{NC mean} + 0.350$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

Q. INTERPRETATION OF RESULTS

Test results are interpreted as ratio of the sample OD450nm and the Cut-Off value (or S/Co) according to the following table:

S/Co	Interpretation
< 0.9	Negative
0.9 - 1.1	Equivocal
> 1.1	Positive

A negative result indicates that the patient has not been infected by HCV or that the blood unit may be transfused. Any patient showing an equivocal result should be tested again on a second sample taken 1-2 weeks later from the patient and examined. The blood unit should not be transfused.

A positive result is indicative of HCV infection and therefore the patient should be treated accordingly or the blood unit should be discarded.

Important notes:

1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgement errors and misinterpretations.
2. Any positive result should be confirmed by an alternative method capable to detect IgG and IgM antibodies (confirmation test) before a diagnosis of viral hepatitis is formulated.
3. As proved in the Performance Evaluation of the product, the assay is able to detect seroconversion to anti HCV core antibodies **earlier** than some other commercial kits. Therefore a positive result, not confirmed with these commercial kits, does not have to be ruled out as a false positive result ! The sample has to be anyway submitted to a confirmation test.
4. As long as the assay is able to detect also IgM antibodies some discrepant results with other commercial products for the detection of anti HCV antibodies - lacking anti hIgM conjugate in the formulation of the enzyme tracer and therefore missing IgM reactivity - may be present. The real positivity of the sample for antibodies to HCV should be then confirmed by examining also IgM reactivity, important for the diagnosis of HCV infection.
5. When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.
6. Diagnosis of viral hepatitis infection has to be done and released to the patient only by a qualified medical doctor.

An example of calculation is reported below:

The following data must not be used instead of real figures obtained by the user.

Negative Control: 0.019 – 0.020 – 0.021 OD450nm
 Mean Value: 0.020 OD450nm Lower than 0.050 – Accepted

Positive Control: 2.189 OD450nm
 Higher than 1.000 – Accepted

Cut-Off = 0.020+0.350 = 0.370
 Calibrator: 0.550 - 0.530 OD450nm
 Mean value: 0.540 OD450nm S/Co = 1.4
 S/Co higher than 1.1 – Accepted

Sample 1: 0.070 OD450nm
 Sample 2: 1.690 OD450nm
 Sample 1 S/Co < 0.9 = negative
 Sample 2 S/Co > 1.1 = positive

R. PERFORMANCES

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications or CTS (art. 5, Chapter 3 of IVD Directive 98/79/EC) and with both incubation procedures (standard and short).

1. LIMIT OF DETECTION

The limit of detection of the assay has been calculated using the short incubation assay by means of the British Working Standard for anti-HCV, NIBSC code 06/188-006. The table below reports the mean OD450nm values of this standard when diluted in negative plasma and then examined.

Dilution	Lot#1	Lot#2	Lot#3
Factor	S/Co	S/Co	S/Co
1 X	3,50	4,00	4,30
2 X	2,10	2,60	2,60
4 X	1,3	1,40	1,30
Negative Plasma	0,25	0,20	0,20

In addition the sample coded Accurun 1 - series 3000 - supplied by Boston Biomedica Inc., USA, has been evaluated "in toto" showing the results below:

Accurun 1 series	Lot#1	Lot#2	Lot#3
Factor	S/Co	S/Co	S/Co
1 X	2,90	3,04	3,40

2. DIAGNOSTIC SPECIFICITY AND SENSITIVITY

The Performance Evaluation of the device was carried out in a trial conducted on more than total 5000 samples.

2.1 Diagnostic specificity

It is defined as the probability of the assay of scoring negative in the absence of specific analyte. A total of more 5000 unselected donors, including 1st time donors, were examined.

The diagnostic specificity was assessed against a kit US FDA approved.

5043 blood donors were tested providing a specificity of 99.5%. 210 hospitalized patients were tested for HCV Ab; a diagnostic specificity of 99.5% was found.

Moreover, diagnostic specificity was assessed by testing 162 potentially interfering specimens (other infectious diseases, E.coli antibody positive, patients affected by non viral hepatic diseases, dialysis patients, pregnant women, hemolized, lipemic, etc.). A value of specificity of 100% was assessed.

No false reactivity due to the method of specimen preparation has been observed. Both plasma, derived with different standard techniques of preparation (citrate, EDTA and heparin), and sera have been used to determine the value of specificity. Frozen specimens have been tested, as well, to check for interferences due to collection and storage. No interference was observed.

2.2 Diagnostic Sensitivity

It defined as the probability of the assay of scoring positive in the presence of specific analyte. The diagnostic sensitivity has been assessed externally on a total number of 348 specimens; a diagnostic sensitivity of 100% was found. Internally more than other 50 positive samples were tested, providing a value of diagnostic sensitivity of again 100%. Positive samples from infections carried out by different genotypes of HCV were tested as well.

Furthermore, most of seroconversion panels available from Boston Biomedica Inc., USA, (PHV) and Zeptometrix, USA, (HCV) have been studied.

Results are reported below for some of them.

Panel	N° samples	Adaltis ¹	Ortho ^{1,2}
PHV 901	11	9	9
PHV 904	7	2	4
PHV 905	9	3	4
PHV 906	7	7	7
PHV 907	7	3	2
PHV 908	13	10	8
PHV 909	3	2	2
PHV 910	5	3	3
PHV 911	5	3	3
PHV 912	3	1	1
PHV 913	4	2	2
PHV 914	9	5	5
PHV 915	4	3	0
PHV 916	8	4	3
PHV 917	10	6	6
PHV 918	8	2	0
PHV 919	7	3	3
PHV 920	10	6	6
HCV 10039	5	2	0
HCV 6212	9	6	7
HCV 10165	9	5	4

Note:

1. Positive samples detected
2. HCV v.3.0

Finally the Product has been tested on the panel EFS Ac HCV, lot n° 06.140817, supplied by the Etablissement Francais Du Sang (EFS), France, with the following results:

EFS Panel Ac HCV

sample	Lot#1 S/Co	Lot#2 S/Co	Lot#3 S/Co	Average value
HCV 1	0,53	0,52	0,55	Negative
HCV 2	3,28	5,91	3,04	Positive
HCV 3	2,17	3,18	2,56	Positive
HCV 4	2,26	2,23	2,35	Positive
HCV 5	6,10	7,06	6,90	Positive
HCV 6	1,66	1,77	1,67	Positive

3. PRECISION

It has been calculated on five samples, one negative and four positives examined in 4 replicates each in six separate runs. The precision, obtained for both protocols, is equivalent.

Results are reported as follows:

Intra-lot results: ElAgen HCV Ab (v.4) Kit - 1st Lot (short incubation assay)

		Precision - %CV		
Sample	S/Co Mean	Within Run	Between Run	Total
Negative	0.03	6.66	10.56	12.48
Positives	1.20	8.52	8.49	12.03
	1.51	7.69	12.22	14.44
	3.57	7.43	11.82	13.97
	11.87	3.42	9.32	9.92

Intra-lot results: ElAgen HCV Ab (v.4) Kit - 1st Lot (long incubation assay)

		Precision - %CV		
Sample	S/Co Mean	Within Run	Between Run	Total
Negative	0.04	4.67	12.34	13.19
Positives	1.47	9.62	11.40	14.92
	1.82	8.92	12.77	15.58
	4.31	4.59	12.88	13.67
	13.78	2.42	8.96	9.26

Inter-lot results: ElAgen HCV Ab (v.4) Kit - 1st, 2nd and 3rd Lot (short incubation assay)

		Precision - %CV		
Sample	Lot 1	Lot 2	Lot 3	
Negative	8,65	8,29	6,13	
Calibrator	4,98	4,44	5,38	
Positive	4,11	3,11	1,37	

The variability shown in the tables above did not result in sample misclassification.

S. SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different ($\pm 50\%$) from OD reported on QC	<ul style="list-style-type: none"> - incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) -incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation) -error in washing or in photometer reading (suggestion: check operating or settings of respective instruments) -contamination of Substrate or Conjugate (suggestion: use only disposable and clean plastic containers)
Low reproducible results	<ul style="list-style-type: none"> -not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes) -error in washing or in reading (suggestion: check operating or settings of respective instruments) -contamination of Substrate (suggestion: use only disposable and clean plastic containers) -pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)
no colorimetric reaction after addition of substrate	<ul style="list-style-type: none"> -some reagent not pipetted - strong contamination of Conjugate or Substrate -errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)
too low reaction (too low ODs)	<ul style="list-style-type: none"> -incubation time too short, incubation temperature too low -incorrect conjugate dilution
too high reaction (too high ODs)	<ul style="list-style-type: none"> -incorrect conjugate dilution -incubation time too long, incubation temperature too high -water quality for wash buffer insufficient (low grade of deionization) -insufficient washing (conjugates not properly removed)
unexplainable outliers	<ul style="list-style-type: none"> -contamination of pipettes, tips or containers -inconstant and insufficient washing (conjugates not properly removed)
too high within-run CV%	<ul style="list-style-type: none"> -reagents and/or strips not pre-warmed to Room Temperature prior to use - plate washer is not washing correctly (suggestion: clean washer head)
too high between-run CV%	<ul style="list-style-type: none"> -incubation conditions not constant (time, temperature) -controls and samples not dispensed at the same time (with the same intervals) (check pipetting order) -person-related variation

T. AUTOMATION

The procedures identified in this Instruction for Use are for manual testing only. When using automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow their approved validation procedures to demonstrate compatibility of this product on automated systems.

U. LIMITATIONS

Repeatable false positive results, not confirmed by RIBA or similar confirmation techniques, were assessed as less than 0.1% of the normal population.

Frozen samples containing fibrin particles or aggregates after thawing have been observed to generate some false results.

BIBLIOGRAPHY

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty MC, Nelson KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB. Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
REFERENCES

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998 ;44 :98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993 ;167 :66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993 ;269 :392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994 ;169 :990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995 ;171 :768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993 ;153 :1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
IRODALOMJEGYZÉK

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. A cikk megtalálható: <http://www.hepnet.com/nih/gretch.html>. A Nemzeti Egészségügyi Hivatal által szervezett, C hepatitisz konferenciára készült, amit 1997 március 24-26.-án tartottak Bethesda Maryland-ban
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
LITERATURA

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, wyd. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turyn, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, wyd. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Waszyngton, DC: US Government Printing Office, 1994; Nr publikacji NIH 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN i in. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH i in. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ i in. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW i in. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB i in. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ i in. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A i in. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK i in. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK i in. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL i in. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW i in. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ i in. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland i A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre i P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. Artykuł dostępny pod adresem: <http://www.hepnet.com/nih/gretch.html>. Napisany w ramach konferencji Narodowego Instytutu Zdrowia dot wirusowego zapalenia wątroby typu C w dniach 24-26 marca 1997 roku w Bethesda w stanie Maryland (USA)
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli i E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams i J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada i T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata i T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human Tlymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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
BIBLIYOGRAFI

1. CDC. Public Health Service inter-agency guidelines for screening donors of blood, plasma, organs, tissues, and semen for evidence of hepatitis B and hepatitis C. *MMWR* 1991;40(No. RR-4):1-17.
2. Alter MJ. Epidemiology of hepatitis C. *Hepatology* 1997;26:62S-5S.
3. McQuillan GM, Alter MJ, Moyer LA, Lambert SB, Margolis HS. A population based serologic study of hepatitis C virus infection in the United States. In Rizzetto M, Purcell RH, Gerin JL, Verme G, eds. *Viral Hepatitis and Liver Disease*, Edizioni Minerva Medica, Turin, 1997, 267-70.
4. Dufour MC. Chronic liver disease and cirrhosis. In Everhart JE, ed. *Digestive diseases in the United States: epidemiology and impact*. US Department of Health and Human Services, Public Health Service, National Institutes of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Washington, DC: US Government Printing Office, 1994; NIH publication no. 94-1447, 615-45.
5. Alter MJ, Hadler SC, Judson FN, et al. Risk factors for acute non-A, non-B hepatitis in the United States and association with hepatitis C virus infection. *JAMA* 1990;264:2231-35.
6. Alter HJ, Holland PV, Purcell RH, et al. Posttransfusion hepatitis after exclusion of commercial and hepatitis-B antigen-positive donors. *Ann Intern Med* 1972;77:691-9.
7. Alter HJ, Purcell RH, Holland PV, Feinstone SM, Morrow AG, Moritsugu Y. Clinical and serological analysis of transfusion-associated hepatitis. *Lancet* 1975;2:838-41.
8. Seeff LB, Wright EC, Zimmerman HJ, McCollum RW, VA Cooperative Studies Group. VA cooperative study of post-transfusion hepatitis and responsible risk factors. *Am J Med Sci* 1975;270:355-62.
9. Feinstone SM, Kapikian AZ, Purcell RH, Alter HJ, Holland PV. Transfusion-associated hepatitis not due to viral hepatitis type A or B. *N Engl J Med* 1975;292:767-70.
10. Choo QL, Kuo G, Weiner AJ, Overby LR, Bradley DW. Isolation of a cDNA clone derived from a blood-borne non-A, non-B viral hepatitis genome. *Science* 1989;244:359-62.
11. Kuo G, Choo QL, Alter HJ, et al. An assay for circulating antibodies to a major etiologic virus of human non-A, non-B hepatitis. *Science* 1989;244:362-4.
12. Alter HJ, Purcell RH, Shih JW, et al. Detection of antibody to hepatitis C virus in prospectively followed transfusion recipients with acute and chronic non-A, non-B hepatitis. *N Engl J Med* 1989;321:1494-1500.
13. Aach RD, Stevens CE, Hollinger FB, et al. Hepatitis C virus infection in post-transfusion hepatitis. An analysis with first- and second-generation assays. *N Engl J Med* 1991;325:1325-9.
14. Alter MJ, Margolis HS, Krawczynski K, Judson, FN, Mares A, Alexander WJ, et al. The natural history of community-acquired hepatitis C in the United States. *N Engl J Med* 1992;327:1899-1905.
15. Alter, MJ. Epidemiology of hepatitis C in the west. *Semin Liver Dis* 1995;15:5-14.
16. Donahue JG, Nelson KE, Muñoz A, et al. Antibody to hepatitis C virus among cardiac surgery patients, homosexual men, and intravenous drug users in Baltimore, Maryland. *Am J Epidemiol* 1991;134:1206-11.
17. Zeldis JB, Jain S, Kuramoto IK, et al. Seroepidemiology of viral infections among intravenous drug users in northern California. *West J Med* 1992;156:30-5.
18. Fingerhood MI, Jasinski DR, Sullivan JT. Prevalence of hepatitis C in a chemically dependent population. *Arch Intern Med* 1993;153:2025-30.
19. Garfein RS, Vlahov D, Galai N, Doherty, MC, Nelson, KE. Viral infections in short-term injection drug users: the prevalence of the hepatitis C, hepatitis B, human immunodeficiency, and human T-lymphotropic viruses. *Am J Pub Health* 1996;86:655-61.
20. Brettler DB, Alter HJ, Deinstag JL, Forsberg AD, Levine PH. Prevalence of hepatitis C virus antibody in a cohort of hemophilia patients. *Blood* 1990;76:254-6.
21. Troisi CL, Hollinger FB, Hoots WK, et al. A multicenter study of viral hepatitis in a United States hemophilic population. *Blood* 1993;81:412-8.
22. Kumar A, Kulkarni R, Murray DL, et al. Serologic markers of viral hepatitis A, B, C, and D in patients with hemophilia. *J Med Virology* 1993;41:205-9.
23. Tokars JI, Miller ER, Alter MJ, Arduino MJ. National surveillance of dialysis associated diseases in the United States, 1995. *ASAIO Journal* 1998;44:98-107.
24. Osmond DH, Charlebois E, Sheppard HW, et al. Comparison of risk factors for hepatitis C and hepatitis B virus infection in homosexual men. *J Infect Dis* 1993;167:66-71.
25. Weinstock HS, Bolan G, Reingold AL, Polish LB: Hepatitis C virus infection among patients attending a clinic for sexually transmitted diseases. *JAMA* 1993;269:392-4.
26. Thomas DL, Cannon RO, Shapiro CN, Hook EW III, Alter MJ. Hepatitis C, hepatitis B, and human immunodeficiency virus infections among non-intravenous drug-using patients attending clinics for sexually transmitted diseases. *J Infect Dis* 1994;169:990-5.
27. Buchbinder SP, Katz MH, Hessel NA, Liu J, O'Malley PM, Alter, MJ. Hepatitis C virus infection in sexually active homosexual men. *J Infect* 1994;29:263-9.
28. Thomas DL, Zenilman JM, Alter HJ, et al. Sexual transmission of hepatitis C virus among patients attending sexually transmitted diseases clinics in Baltimore--an analysis of 309 sex partnerships. *J Infect Dis* 1995;171:768-75.
29. Thomas DL, Factor SH, Kelen GD, Washington AS, Taylor E Jr, Quinn TC. Viral hepatitis in health care personnel at The Johns Hopkins Hospital. *Arch Intern Med* 1993;153:1705-12.
30. Cooper BW, Krusell A, Tilton RC, Goodwin R, Levitz RE. Seroprevalence of antibodies to hepatitis C virus in high-risk hospital personnel. *Infect Control Hosp Epidemiol* 1992;13:82-5.
31. Abdel-Hamid, M., M. El-Day, S. El-Kafrawy, N. Mikhail, G.T. Strickland, and A.D. Fix. 2002. Comparison of second- and third-generation enzyme immunoassays for detecting antibodies to hepatitis C virus. *J. Clin. Microbiol.* 40:1656-1659.
32. Dusheiko, G., H. Schlimovitz-Weiss, D. Brown, F. McOmish, P.-L. Yap, S. Sherlock, N. McIntyre, and P. Simmonds. 1994. Hepatitis C virus genotypes: an investigation of type-specific differences in geographic origin and disease. *Hepatology* 19:13-18.
33. Gretch, D. Diagnostic tests for hepatitis C. The article can be found at: <http://www.hepnet.com/nih/gretch.html>. It was written as part of a National Institute of Health Conference on Hepatitis C, held from March 24-26, 1997 in Bethesda Maryland
34. Mondelli, M.U., A. Cerino, F. Bono, A. Cividini, A. Maccabruni, M. Aricò, A. Malfitano, G. Barbarini, V. Piazza, L. Minoli, and E. Silini. 1994. Hepatitis C virus (HCV) core serotype in chronic HCV infection. *J. Clin. Microbiol.* 32:2523-2527.
35. Ohno, T., M. Mizokami, R.-R. Wu, M.G. Saleh, K.-I. Ohba, E. Orito, M. Mukaide, R. Williams, and J.Y.N. Lau. 1997. New hepatitis C virus (HCV) genotyping system that allows for the identification of HCV genotype 1a, 1b, 2a, 2b, 3a, 3b, 4, 5a, and 6a. *J. Clin. Microbiol.* 35:201-207.
36. Takada, N., S. Takase, N. Enomoto, A. Takada, and T. Date. 1992. Clinical backgrounds of the patients having different types of hepatitis C virus genome. *J. Hepatol.* 14:35-40.
37. Yoshioka, K., S. Kakumu, T. Wakita, T. Ishikawa, Y. Itoh, M. Takayanagi, Y. Higashi, M. Shibata, and T. Morishima. 1992. Detection of hepatitis C virus by polymerase chain reaction and response to interferon- α therapy: relationships to genotypes of hepatitis C virus. *Hepatology* 16:293-299.

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**A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUANTITATIVE DETERMINATION
OF TOTAL IgG IN HUMAN BIOLOGICAL FLUIDS**

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
Authorized Representative in EU:
Polmed.de
Steinacker 20, D-73773
Aichwald, Germany
e-mail: info@polmed.de

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* Blood serum or plasma

For other tested materials, see table M.

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CONTENT

1. INTENDED USE	10
2. SUMMARY AND EXPLANATION	10
3. PRINCIPLE OF THE TEST	10
4. WARNINGS AND PRECAUTIONS	11
5. KIT COMPONENTS	12
6. SPECIMEN COLLECTION AND STORAGE	13
7. TEST PROCEDURE	13
8. QUALITY CONTROL	15
9. CALCULATION OF RESULTS	15
10. EXPECTED VALUES	16
11. PERFORMANCE CHARACTERISTICS	16
12. LITERATURE	16

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IgG IN HUMAN BIOLOGICAL FLUIDS

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of total IgG in biological fluids.

This kit is designed for measurement of total IgG in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Immunoglobulin G (IgG) is the main part of serum γ - globulin fraction. IgG is secreted during secondary immune response and plays a key role in humoral immunity. Decrease of serum IgG concentration below 5 g/l is a marker of severe life-threatening immunodeficiency. Determination of serum IgG concentration and IgG/IgA/IgM ratios can be used for monitoring of humoral immune status. Marked elevation of serum IgG may be observed in chronic inflammation, autoimmune diseases and myeloma.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgG-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies - murine monoclonal to human total IgG, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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4.4. Avoid contact with stop solution containing 5.0% H_2SO_4 . It may cause skin irritation and burns.

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5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP total IgG EIA strips, 8x12 wells	1	pcs		until exp.date
2	CAL 1-5 Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 1; 5; 10; 25 g/l	5	pcs	blue (Cl - colourless)	2 months
3	CONTROL Control serum (1 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 14 ml	1	pcs	red	until exp.date
5	DIL SPE EIA sample buffer 100 ml	1	pcs	blue	until exp.date
6	SUBS TMB Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	BUF WASH 26X Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate - until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8	STOP Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	N003 Plate sealing tape	2	pcs		N/A
10	K2711 Instruction total IgG EIA	1	pcs		N/A
11	K271Q QC data sheet total IgG EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 10–250 µl;
- Dry thermostat for 37 °C ±0.1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

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Specimens may be stored for up to 48 hours at +2...+8 °C before testing.

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- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

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7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA sample buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA sample buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1-5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the unknown diluted sample (DILUTION 2) into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5-7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at +37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape.
8	Incubate 30 minutes at +37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10-20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA sample buffer into the well, µl	Sample into the well, µl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 µl of sample + 990 µl of diluent = DILUTION 1. 10 µl of DILUTION 1 + 490 µl of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		90	10	0.002
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 µl of sample + 500 µl of diluent	0	100	0.01

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
9.2. Plot a calibration curve on graph paper: OD versus total IgG concentration.

9.3. Determine the corresponding concentration of total IgG in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is

recommended due to non-linear shape of curve.

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10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgG. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
newborn	7.0	15
1-3 month	2.7	8.0
4-6 month	1.8	8.5
7-12 month	3.5	12
1-6 yrs	6.5	18
7-11 yrs	8.5	15
> 11 yrs	9.0	20

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

11.3. Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgG concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgG concentrations. The recovery percentages ranged from 90 to 110%.

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Russian Diagnostic
Manufacturers Association

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Russian Association
of Medical Laboratory
Diagnosticians

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IGA IN HUMAN BIOLOGICAL FLUIDS

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of total IgA in biological fluids.

This kit is designed for measurement of total IgA in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Immunoglobulin A (IgA) is a main factor of mucosal immune response to bacteria and viruses. Selective IgA deficiency is one of the most frequent hereditary disorders causing chronic infections inflammation in gastrointestinal, urinary and respiratory systems. Determination of IgA concentration in serum and other biological fluids can be used as screening for selective IgA deficiency and other immunodeficiency syndromes. Marked elevation of serum IgA is observed in some autoimmune diseases and IgA myeloma.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgA-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human total IgA, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5.0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	SORB MTP total IgA EIA strips, 8x12 wells	1	pcs		until exp.date
2	CAL 1-5 Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 0.1; 0.5; 2; 5 g/l	5	pcs	blue (C1 – colourless)	2 months
3	CONTROL Control serum (1 ml)	1	pcs	colourless	2 months
4	CONJ HRP Conjugate, 14 ml	1	pcs	blue	until exp.date
5	DIL SPE EIA buffer 100 ml	1	pcs	blue	until exp.date
6	SUBS TMB Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	BUF WASH 26X Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8	STOP Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	N003 Plate sealing tape	2	pcs		N/A
10	K275I Instruction total IgA EIA	1	pcs		N/A
11	K275Q QC data sheet total IgA EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 μl , is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250 μl ;
- Dry thermostat for $+37\text{ }^{\circ}\text{C} \pm 0.1\text{ }^{\circ}\text{C}$
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at $+2\dots+8\text{ }^{\circ}\text{C}$ upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at $+2\dots+8\text{ }^{\circ}\text{C}$ before testing.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature ($+18\dots+25\text{ }^{\circ}\text{C}$) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1-5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1-5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the diluted unknown sample into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5-7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at +37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells.
8	Incubate 30 minutes at +37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10-20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA buffer into the well, μl	Sample into the well, μl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 μl of sample + 990 μl of diluent = DILUTION 1. 10 μl of DILUTION 1 + 490 μl of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	5 μl of sample + 500 μl of diluent	90	10	0.2
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		80	20	0.001
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus total IgA concentration.

9.3. Determine the corresponding concentration of total IgA in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 g/l	0.08
CAL 2	0.1 g/l	0.40
CAL 3	0.5 g/l	1.32
CAL 4	2 g/l	2.45
CAL 5	5 g/l	2.92

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgA. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
Healthy donors	0.9	5.0
>61 yr	1.0	6.5
newborn	-	0.05
1-3 month	0.06	0.6
4-6 month	0.1	1.0
7-12 month	0.35	1.7
1-6 yrs	0.8	2.2
7-11 yrs	0.9	2.6

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
IgG	<0.1
IgM	<0.1
IgE	<0.1

11.2. Analytical sensitivity. Sensitivity of the assay was assessed as being 0.06 g/l.



11.3. Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgA concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgA concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

1. Heiddis B. Valdimarsdottir and Arthur A. Stone – Psychosocial Factors and Secretory Immunoglobulin A. Critical Reviews in Oral Biology & Medicine, Jan 1997; 8: 461 – 474.

2. Amir H Abdul Latiff and Michael A Kerr – The clinical significance of immunoglobulin A deficiency. Ann Clin Biochem, Mar 2007; 44: 131 – 139.

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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF TOTAL IgM IN HUMAN BIOLOGICAL FLUIDS

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of total IgM in biological fluids.

This kit is designed for measurement of total IgM in biological fluids. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 42 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Immunoglobulin M (IgM) is secreted during primary immune response and exists in monomeric and pentameric forms. Elevated serum IgM is observed in chronic inflammation, macroglobulinemia and IgM myeloma. Decreased IgM level may occur in some immunodeficiency syndromes.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to human total IgM-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Unbound material is removed by washing procedure. Second antibodies – murine monoclonal to human total IgM, labelled with peroxidase enzyme, are then added into the microwells. After subsequent washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

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5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	total IgM EIA strips, 8x12 wells	1	pcs		until exp.date
2	total IgM EIA strips, 8x12 wells Calibrator set, 1 ml each. The set contains 5 calibrators: 0; 0.5; 2; 5; 10 g/l	5	pcs	purple (C1 – colourless)	2 months
3	Control serum (1 ml)	1	pcs	colourless	2 months
4	Conjugate, 14 ml	1	pcs	purple	until exp. date
5	EIA buffer 100 ml	1	pcs	blue	until exp. date
6	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate – until exp.date Diluted washing solution – 45 days at 2-8 °C or 15 days at RT
8	Stop solution, 14 ml	1	pcs	colourless	until exp. date
9	Plate sealing tape	2	pcs		N/A
10	Instruction total IgM EIA	1	pcs		N/A
11	QC data sheet total IgM EIA	1	pcs		N/A

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- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 5–250 µl;
- Dry thermostat for 37 °C ±0.1 °C
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

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After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

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7.4. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1–5 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	Dilute samples using buffer DIL SPE (EIA sample buffer) 5000 fold. See table M for dilution modes and factors for different types of analyzed material. Do not dilute control sample and calibrators.
3	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA sample buffer). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
4	Pipet 100 µl of calibrators CAL 1–5 and control samples CONTROL into allocated wells. For testing of blood serum or plasma pipet 100 µl of the unknown diluted sample (DILUTION 2) into the allocated wells. See table M for the volumes of other materials. Pipetting should be made within 3 minutes, to ensure an uniform incubation time for all samples. Carefully mix the contents of the wells by short horizontal rotating of the plate for 5–7 seconds and cover the wells by plate adhesive tape (included into the kit).
5	Incubate 30 minutes at 37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate BUF WASH 26X by distilled water. Minimal quantity of washing solution should be 250 µl per well. Wash strips 3 times.
7	Dispense 100 µl of CONJ HRP into the wells.
8	Incubate 30 minutes at 37 °C.
9	Wash the strips 5 times.
10	Dispense 100 µl of SUBS TMB into the wells.
11	Incubate 10–20 minutes at +18...+25 °C.
12	Dispense 100 µl of STOP into the wells.
13	Measure OD (optical density) at 450 nm.
14	Set photometer blank on first calibrator.
15	Apply point-by-point method for data reduction. Use Calculation factor listed in table M to calculate analyte concentration in different material types.

7.5. Sample processing

Material type	Notes on material collection, storage and handling	Sample dilution example	EIA sample buffer into the well, µl	Sample into the well, µl	Calculation factor
blood serum or plasma	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.	10 µl of sample + 990 µl of diluent = DILUTION 1. 10 µl of DILUTION 1 + 490 µl of diluent = DILUTION 2	0	100	1
saliva	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		90	10	0.002
urine	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		50	50	0.0004
cerebrospinal fluid	Grossly hemolytic, lipemic, or turbid samples should be avoided and should be treated by centrifugation before testing.		80	20	0.001

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
9.2. Plot a calibration curve on graph paper: OD versus total IgM concentration.

9.3. Determine the corresponding concentration of total IgM in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

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10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgM. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, g/l	
	Lower limit	Upper limit
Healthy donors	0.7	3.7
newborn	0.1	0.35
1-3 month	0.12	0.9
4-6 month	0.25	1.2
7-12 month	0.35	1.0
1-6 yrs	0.55	2.2
7-11 yrs	0.65	1.7

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11.3. Linearity. Linearity was checked by assaying dilution series of 5 samples with different total IgM concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery. Recovery was estimated by assaying 5 mixed samples with known total IgM concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

Erik J. Wiersma, Cathy Collins, Shafie Fazel, and Marc J. Shulman Structural and Functional Analysis of J Chain-Deficient IgM J. Immunol., Jun 1998; 160: 5979 – 5989.

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xemahelp

xemahelp@gmail.com

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
total IgE in human serum or plasma

Total IgE EIA

Catalogue number **REF** **K200**

For 96 determinations

In vitro diagnostic medical device


XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.com.ua

Authorized Representative in EU:
Polmed.de Beata Rozwadowska
Fichtenstr. 12A, 90763 Fuerth, Germany
tel.:+ 49 911 931 639 67
E-mail: info@polmed.de
www.polmed.de

ASSAY PROCEDURE

K200

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	8
12. EXPECTED VALUES	8
13. PERFORMANCE CHARACTERISTICS	9
14. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
total IgE in human serum or plasma
Total IgE EIA

1. INTENDED USE

The Total IgE EIA kit is an enzyme immunoassay, intended for the quantitative determination of total IgE concentration in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Total immunoglobulin E (IgE) serum level is widely reported as the laboratory marker of atopic diseases such as atopic asthma, atopic dermatitis, and pollenosis. An atopic (IgE-dependent) mechanism can also underlie gastroenterocolitis, urticaria, other forms of vasculitis (including systemic), cholecystitis, vulvovaginitis, and cystitis. Part of the drug allergy (mainly to penicillin and protein drugs) also develops according to the IgE-dependent mechanism. In all of the conditions listed above, the production of high titers of specific IgE antibodies can lead to an increase in the level of total IgE in the serum. A particularly high level of total IgE is characteristic of atopic dermatitis. In addition to atopic diseases, total serum IgE is significantly increased in parasitic infestations and mycoses (especially systemic), rarely in systemic autoimmune diseases and immunodeficiency states (especially in hyper-IgE syndrome), as well as in mastocytosis (mast cell tumor) and extremely rare IgE-myeloma. A decrease in the level of total IgE in serum (below 15 IU/ml in adults) is a rare and little-studied phenomenon described in hypogammaglobulinemia, some autoimmune diseases, ulcerative colitis, and primary biliary cirrhosis.

3. TEST PRINCIPLE

The determination of the total IgE is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to human IgE. Second antibodies – rabbit polyclonal antibodies to IgE conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes three stages of incubation:

- during the first stage the total IgE from the specimen is captured by the monoclonal antibodies coated onto the microwell surface;
- during the second stage horseradish peroxidase-conjugated with rabbit polyclonal antibodies bind to free epitopes of immobilized total IgE, fixed in the formed at the previous stage complexes;
- during the third stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured total IgE in the serum specimen (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of total IgE in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P200Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to total IgE; ready to use
C200Z	CAL 1	Calibrator C1	0.8 mL	1	Solution based on phosphate buffer, free of total IgE, with preservative, ready to use (yellow liquid)
C200Z	CAL 1-5	Calibrators	0.8 mL	5	Solutions based on phosphate buffer, containing 50; 200; 500 and 1000 IU/mL of total IgE, ready to use (red liquids)
Q200Z	CONTROL	Control serum	0.8 mL	1	Solution based on human serum, containing of known total IgE content, with preservative, ready to use (colourless liquid)
T200Z	CONJ HRP	Conjugate Solution	14 mL	1	Solution of rabbit polyclonal antibodies to human total IgE conjugated to the horseradish peroxidase; ready to use (red liquid)
S011Z	DIL	EIA Buffer	14 mL	1	Buffer solution with detergent and preservative, ready to use (blue liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	30 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (3 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength or 450\620-680 nm;
- dry thermostat for 37 °C±2 °C;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below; do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The Total IgE EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The Total IgE EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- EIA Buffer, Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing solution preparation

Add the contents of the 30 mL washing solution concentrate vial to 750 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of washing solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the washing solution concentrate, mL	2.5	5	7.5	10	12.5	15	17.5	20	22.5	25	27.5	30
Volume of water, mL	62.5	125	187.5	250	312.5	375	437.5	500	562.5	625	687.5	750

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-5) and 2 wells for Control Serum (Q)).
- 10.2 Dispense **50 µL of EIA Buffer** to all wells.
- 10.3 Dispense **50 µL of Calibrators and Control Serum as well as 50 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators, and Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5	SAMP13	SAMP13						
D	CAL4	CAL4	SAMP6	SAMP6	SAMP14	SAMP14						
E	CAL5	CAL5	SAMP7	SAMP7	SAMP15	SAMP15						
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.4 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **30 minutes at +37 °C**.
- 10.5 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 3 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.6 Add **100 µL of Conjugate Solution** to all wells.
- 10.7 Cover strips with a plate sealing tape and incubate for **30 minutes at +37 °C**.
- 10.8 At the end of the incubation period, aspirate and wash each well 5 times as described in 10.5.
- 10.9 Add **100 µL of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 20 minutes**.
- 10.10 Add **100 µL of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.11 Read the optical density (OD) of the wells at 450nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the stop solution. Set photometer blank on CAL1.
- 10.12 Plot a calibration curve in linear coordinates: (x) is the concentration of total IgE in the Calibrators IU/mL, (y) – OD versus concentration of total IgE (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.
- 10.13 Determine the corresponding concentration of total IgE in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

12.1. Therapeutical consequences should not be based on the results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for total IgE. Based on data obtained by XEMA LLC, the following normal range is recommended (see below).

NOTE: values of total IgE concentrations in the tested samples that are below the LoD (3 IU/mL) and also exceed the value of the upper calibrator (1000 IU/mL) should be provided in the following form : «the total IgE concentration of tested sample X is «lower than 3 IU/mL» or «higher than 1000 IU/mL».

12.2. The calibrators concentration values of the Total IgE EIA kit are expressed in IU/mL. To calculate concentrations in ng/mL, the received concentration value in IU/mL shall be multiplied by 2.4.

$$1 \text{ IU/mL} = 2.4 \text{ ng/mL.}$$

Sex, age	Units, IU/mL		Units alternative, ng/mL	
	Lower limit	Upper limit	Lower limit	Upper limit
< 6 months	-	12	-	28.8
6-12 months	-	30	-	72.0
1-3 yrs	-	45	-	108.0
4-6 yrs	-	70	-	168.0
7-9 yrs	-	90	-	216.0
10-15 yrs	-	120	-	288.0
>15 yrs	-	130	-	312.0

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, IU/mL	CV, %
1	10.6	4.33
2	116.2	5.47

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, IU/mL	CV, %
1	12.5	8.36
2	113.4	1.47

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, IU/mL	Concentration2, IU/mL	Concentration3, IU/mL	CV, %
1	12.7	13.3	12.3	3.66
2	115.5	117.8	115.1	1.25

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known total IgE concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 50-1000 IU/mL $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest total IgE concentration in the serum or plasma sample that is detected by the Total IgE EIA kit is no lower than 3 IU/mL.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for Total IgE EIA kit is 50IU/mL.

13.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of total IgE with other analytes is shown in the table:

Analyte	Concentration, IU/mL	Cross-reactivity, %
IgA	1000	Not detected
IgM	1000	Not detected
IgG	1000	Not detected

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
thyroid stimulating hormone
in human serum or plasma

TSH EIA

Catalogue number **REF** **K201**

For 96 determinations

In vitro diagnostic medical device


XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Authorized Representative in EU:
Polmed.de Beata Rozwadowska
Fichtenstr. 12A, 90763 Fuerth, Germany
tel.:+ 49 911 931 639 67
E-mail: info@polmed.de
www.polmed.de


ASSAY PROCEDURE

Conjugate Solution

100 μ L


Dispensing of calibrators, control serum and test samples

50 μ L

SORB MTP

Incubation 1


+18...+25 °C

90'

Washing 5 times


OD measuring, calculation of results


450/620-680 nm

Stop Solution

100 μ L


Incubation 2

+18...+25 °C

20'

Substrate Solution

100 μ L

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	3
4. KIT COMPONENTS	4
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	5
6. WARNING AND PRECAUTIONS	5
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	6
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	6
9. REAGENTS PREPARATION	7
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	9
12. EXPECTED VALUES	9
13. PERFORMANCE CHARACTERISTICS	9
14. REFERENCES	11
SAMPLES IDENTIFICATION PLAN	12

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
thyroid stimulating hormone
in human serum or plasma
TSH EIA

1. INTENDED USE

The TSH EIA kit is an enzyme immunoassay, intended for the quantitative determination of thyroid stimulating hormone in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Thyroid stimulating hormone (TSH) is a glycoprotein with molecular weight ca.30 kDa which is secreted by hypophysis. A molecule of TSH consists of two noncovalently bound subunits: α and β . β -subunit determines biological activity and immunological specificity of TSH.

TSH stimulates thyroid gland to secrete thyroid hormones. When the concentration of these hormones in blood serum increases secretion of TSH is inhibited; on the contrary, when the level of thyroid hormones decreases, in the pituitary gland, the release of TSH increases, and therefore the production and release increases thyroid hormones. TSH secretion is subject to circadian rhythms with highest levels seen early in the morning (6 a.m.). Changes of TSH blood level during a day are not significant; nevertheless, if the results do not correspond with clinical status and other laboratory data, it is recommended to take and test another blood sample.

Determination of TSH level in serum is recommended in the following states and conditions:

- 1) diagnostics of dysfunction of the thyroid gland;
- 2) hypothyroidism (TSH level is increased. The diagnosis is confirmed by low concentrations of total and free T4 and T3. In mild subclinical forms when T4 and T3 levels are within normal ranges, determination of TSH concentration is critical);
- 3) hyperthyroidism (synthesis and secretion of TSH are inhibited); monitoring of replacement therapy;
- 4) screening for congenital hypothyroidism (on the fifth day of life, the level is determined TSH in a blood spot on filter paper or in blood serum). TSH level elevated at birth (up to 35 mIU/L), but after a few days it decreases to basal (both in boys and in girls).

Serum TSH level is elevated during pregnancy, after physical stress, in individuals with lowered blood pressure and lowered temperature. Secretion of TSH is inhibited by Cortisol and Growth hormone. Low TSH levels are often seen in elderly people, in patients with chronic renal insufficiency, liver cirrhosis, in retardation of sexual development, in secondary amenorrhea, Cushing syndrome, acromegaly.

3. PRINCIPLE OF THE TEST

The determination of TSH is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to β -chain of human TSH. Second antibodies – Fab 2 fragment of murine monoclonal antibodies to human TSH conjugated to the horseradish peroxidase is used as enzyme conjugate.

The analysis procedure includes two stages of incubation:

- during the first stage TSH from the specimen is captured by the antibodies coated onto the microwell surface, as well as horseradish peroxidase-conjugated monoclonal antibodies bind to free epitopes of immobilized TSH;

- during the second stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured TSH in the serum specimen (plasma). The concentration is determined according to the calibration graph of the dependence of the optical density on the content of TSH in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P201Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to β -chain of human TSH; ready to use
C201Z	CAL 1	Calibrator C1	2 mL	1	Solution based on phosphate buffer (pH 7.2-7.4), free of human TSH, with preservative, ready to use (yellow liquid)
C201Z	CAL 2-6	Calibrators	0.8 mL	5	Solution based on phosphate buffer (pH 7.2-7.4), containing 0,2; 1; 5; 10 and 20 mIU/L of human TSH, with preservative, ready to use (red liquids)
Q201Z	CONTROL	Control Serum	0.8 mL	1	Solution based on human serum, containing of known human TSH content, with preservative, ready to use (colourless liquid)
T201Z	CONJ HRP	Conjugate Solution	14 mL	1	Solution of Fab 2 fragment of murine monoclonal antibodies to human TSH conjugated to the horseradish peroxidase; ready to use (blue liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)

The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength or 450\620-680 nm;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The TSH EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The TSH EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months.

NOTE: Single freezing of Calibrators and Control Serum in aliquots is allowed.

- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of Washing Solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

9.4. Samples preparation

If suggested analyte concentration in the sample exceeds the 20 mIU/L, additionally dilute this sample accordingly, using (Calibrator C1). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.

NOTE: in order to obtain reliable results, we recommend to use several successive dilutions of the blood serum (plasma) sample

Do not dilute Control Serum and Calibrators!

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- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 14 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-6) and 2 wells for Control Serum (Q)).
- 10.2 If necessary, dilute the test samples as described in 9.4.
- 10.3 Dispense **100 µL of Conjugate Solution** to all wells.
- 10.4 Dispense **50 µL of Calibrators and Control Serum as well as 50 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators, and Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP2	SAMP2	SAMP10	SAMP10						
B	CAL2	CAL2	SAMP3	SAMP3	SAMP11	SAMP11						
C	CAL3	CAL3	SAMP4	SAMP4	SAMP12	SAMP12						
D	CAL4	CAL4	SAMP5	SAMP5								
E	CAL5	CAL5	SAMP6	SAMP6								
F	CAL6	CAL6	SAMP7	SAMP7								
G	Q	Q	SAMP8	SAMP8								
H	SAMP1	SAMP1	SAMP9	SAMP9								

- 10.5 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **90 minutes at room temperature (+18...+25°C)**.
- 10.6 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the wash solution volume can be increased to 350 µL.
- 10.7 Add **100 µL of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 20 minutes**.
- 10.8 Add **100 µL of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.9 Read the optical density (OD) of the wells at 450nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the Stop Solution. Set photometer blank on CAL1.
- 10.10 Plot a calibration curve in linear coordinates: (x) is the TSH concentration in the calibrators mIU/L, (y) – OD versus TSH concentration (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.
- 10.11 Determine the corresponding concentration of TSH in tested samples from the calibration curve. In the case of preliminary dilution of the test sample (see 9.4), the obtained result should be multiplied by the dilution factor.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on the results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for TSH. Based on data obtained by XEMA, the following normal range is recommended (see below).

NOTE: values of TSH concentrations in the tested samples that are below the LoD (0.04 mIU/L) and also exceed the value of the upper Calibrator (20 mIU/L) should be provided in the following form : «the TSH concentration of tested sample X is «lower than 0.04 mIU/L» or «higher than 20 mIU/L».

Sex, age	Units, mIU/L	
	Lower limit	Upper limit
Healthy donors	0.3	4.0

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, mIU/L	CV, %
1	2.12	7.2
2	3.64	3.8

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, mIU/L	CV, %
1	2.27	12.0
2	3.87	6.4

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, mIU/L	Concentration2, mIU/L	Concentration3, mIU/L	CV, %
1	2.32	2.02	1.81	9.9
2	3.71	3.56	3.32	5.6

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known TSH concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 0.2-10 mIU/L \pm 10%.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest TSH concentration in the serum or plasma sample that is detected by the TSH EIA kit is no lower than 0.04 mIU/L.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for TSH EIA kit is 0,15 mIU/L.

13.1.5 Hook Effect

Hook effect is absent for all samples up to reasonably foreseen concentrations 20 mIU/L.

13.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of TSH with other analytes is shown in the table:

Analyte	Cross-reactivity, %
HCG	< 0.1
LH	< 0.1
FSH	< 0.1

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
triiodothyronine in human serum or plasma

T3 EIA

Catalogue number **REF** **K211**

For 96 determinations

In vitro diagnostic medical device

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.com.ua

Authorized Representative in EU:
Polmed.de Beata Rozwadowska
Fichtenstr. 12A, 90763 Fuerth, Germany
tel.:+ 49 911 931 639 67
E-mail: info@polmed.de
www.polmed.de


ASSAY PROCEDURE

Preparation of working conjugate solution

Conjugate Dilution Buffer + Conjugate Concentrate = Working conjugate solution


1 : 1

450 μ L 450 μ L 450 μ L


Dispensing of Calibrators, Control Serum and test samples

test samples 25 μ L SORB MTP

Working Conjugate Solution

100 μ L


Incubation 1

+37 $^{\circ}$ C


60'

continuous shaking


OD measuring, calculation of results

450 nm

Stop Solution


100 μ L

Incubation 2

+18-25 $^{\circ}$ C


15'


Substrate Solution

100 μ L

SUBS TMB

Washing 5 times

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	6
11. TEST VALIDITY	7
12. EXPECTED VALUES	8
13. PERFORMANCE CHARACTERISTICS	8
14. REFERENCES	9
SAMPLES IDENTIFICATION PLAN	10

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
triiodothyronine in human serum or plasma
T3 EIA

1. INTENDED USE

The T3 EIA kit is an enzyme immunoassay, intended for the quantitative determination of triiodothyronine in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Triiodothyronine (T3) is a hormone with a molecular weight of 651 Da, 58% of which is iodine. Thyroid hormones thyroxin (T4) and 3,5,3'-triiodothyronine (T3) exert regulatory influences on growth, differentiation, cellular metabolism and development of skeletal and organ systems. T4 and T3 in blood are found both in free and bound form – mostly, they are bound to thyroxin binding globulin (TBG). Only free forms of T3 and T4 exert hormonal activity also their percentage is very low – 0.3% for T3 and 0.03% for T4.

The concentration of T3 is much less than that of T4 but its metabolic activity is about 3 times greater. About 80% of T3 is produced in peripheral tissues by deiodination of T4, and only 20% is secreted by thyroid gland. That is why in hypothyroid patients T3 level may for a long time remain on the lower limit of the normal range, because its loss may be compensated by enhanced conversion of T4 into T3.

Determination of T3 level is most useful in T3-hyperthyroidism because 5-10% of such patients do not show significant changes in T4 level while concentration of T3 is highly elevated. Elevated T3 levels are seen in early thyroid hypofunction, after intake of estrogens, oral contraceptives, heroin, methadone, during pregnancy.

Decreased concentrations of T3 are found in initial stage of hyperthyroidism, acute and subacute thyroiditis, after intake of androgens, dexamethasone, salicylates. Decreased concentrations of T3 are found in initial stage of hyperthyroidism, acute and subacute thyroiditis, after intake of androgens, dexamethasone, salicylates.

3. TEST PRINCIPLE

The determination of triiodothyronine is based on the competition principle of the enzyme immunoassay. On the inner surface of the microplate wells are immobilized specific rabbit polyclonal to T3 antibodies. T3 conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage T3 from the specimen competes with the conjugated T3 for coating antibodies. As a result, a complex bounded to the solid phase and containing peroxidase is formed.

- during the second stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is inversely related to the quantity of the measured T3 in the serum specimen (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of T3 in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P211Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with rabbit polyclonal antibodies to T3, ready to use;
C211Z	CAL 1	Calibrator C1	0.5 mL	1	Solution based on tris buffer (pH 7.2-7.4), free of T3, with preservative, ready to use (yellow liquid)
C211Z	CAL 2-5	Calibrators	0.5 mL	4	Solutions based on tris buffer (pH 7.2-7.4), containing 0,75; 1,5; 7,5 and 15 nmol/L of T3, with preservative, ready to use (blue liquids)
Q211Z	CONTROL	Control serum	0.5 mL	1	Solution based on human plasma, containing of known T3 content, with preservative, ready to use (colourless liquid)
T211XZ	CONJ 2X	Conjugate Concentrate	7 mL	1	Solution of T3 conjugated to the horseradish peroxidase; 2x concentrate (purple liquid)
ST211Z	DIL CONJ	Conjugate Dilution Buffer	7 mL	1	Buffer solution with detergent ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	14 ml (MЛ)	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 ml (MЛ)	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 ml (MЛ)	1	5,0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength;
- shaker maintaining a speed of 500 rpm for +37 °C±2°C;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below; do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The T3 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The T3 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Concentrate, Conjugate Dilution Buffer, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of Washing Solution concentrate and dilute it 26 times with distilled or deionized water.

9.4. Working conjugate solution preparation

Prepare a working conjugate solution by 2 dilutions of Conjugate Concentrate in Conjugate Dilution Buffer (eg, 450 µL of concentrate + 450 µL of Conjugate Dilution Buffer). In the case of partial use of the kit, take the necessary amount of Conjugate Concentrate and dilute it 2 times with Conjugate Dilution Buffer, since the working conjugate solution in a diluted form is not stored for a long time.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550
Volume of Conjugate Concentrate, mL	0.45	0.9	1.35	1.8	2.25	2.7	3.15	3.6	4.05	4.5	4.95	5.4
Volume of Conjugate Dilution Buffer, mL	0.45	0.9	1.35	1.8	2.25	2.7	3.15	3.6	4.05	4.5	4.95	5.4

10. ASSAY PROCEDURE

- Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each calibrator (CAL 1-5) and 2 wells for control serum (Q)).
- Prepare Working conjugate solution as described in 9.4.
- Dispense **25 µL of Calibrators and Control Serum as well as 25 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

Note: during performing several independent series of tests, Calibrators, and Control Sample should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5								
D	CAL4	CAL4	SAMP6	SAMP6								
E	CAL5	CAL5	SAMP7	SAMP7								
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.4 Dispense **100 µL of Working conjugate solution** to all wells.
- 10.5 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C with continuous shaking 500 rpm**.
- 10.6 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.7 Add **100 µL of Substrate Solution** to all wells. The introduction of the substrate solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.8 Add **100 µL of Stop Solution** to all wells in the same order as the substrate solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.9 Read the optical density (OD) of the wells at 450nm using a microplate photometer within 5 minutes of adding the Stop Solution.
- 10.10 Plot a calibration curve in semi-logarithmic coordinates: (x) is the decimal logarithm of the T3 concentration in the calibrators nmol/L, (y) – OD versus T3 concentration (OD 450 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve. Adjust the concentration of CAL1 to an infinitesimally small value, for example, 0.001 nmol/L.
- 10.11 Determine the corresponding concentration of T3 in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 1.2, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for T3. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE: values of T3 concentrations in the tested samples that are below the LoD (0.2 nmol/L) and also exceed the value of the upper calibrator (15 nmol/L) should be provided in the following form: «the T3 concentration of tested sample X is «lower than 0.2 nmol/L» or «higher than 15 nmol/L».

The concentration values of the T3 EIA kit calibrators are expressed in nmol/L. To convert the concentration in ng/mL it is necessary to multiply by 0.65 the obtained concentration value in nmol/L.

$$1 \text{ nmol/L} = 0.65 \text{ ng/mL}$$

Sex, age	Units, nmol/L		Units alternative, ng/mL	
	Lower limit	Upper limit	Lower limit	Upper limit
Healthy donors	1.2	3.2	0.8	2.1

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, nmol/L	CV, %
1	2.32	9.16
2	1.45	9.66

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, nmol/L	CV, %
1	1.38	9.89
2	1.75	8.41

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, nmol/L	Concentration2, nmol/L	Concentration3, nmol/L	CV, %
1	2.12	2.02	2.27	13.9
2	1.56	1.44	1.81	15.6

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known T3 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 0.75 –15 nmol/L $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest T3 concentration in the serum or plasma sample that is detected by the T3 EIA kit is no lower than 0.2 nmol/L.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for T3 EIA kit is 0.55 nmol/L.

3.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of T3 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
L-Thyroxin	0.01
D-Thyroxin	0.04

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.com.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
thyroxin in human serum or plasma

T4 EIA

Catalogue number **REF** **K212**

For 96 determinations

In vitro diagnostic medical device

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Authorized Representative in EU:
Polmed.de Beata Rozwadowska
Fichtenstr. 12A, 90763 Fuerth, Germany
tel.:+ 49 911 931 639 67
E-mail: info@polmed.de
www.polmed.de

ASSAY PROCEDURE

Dispensing of Calibrators, Control Serum and test samples

test samples
CAL
CONTROL
25 µL
SORB MTP

Conjugate Solution

CONJ HRP
100 µL

Incubation 1

+37 °C
60'

Washing 5 times

Washing 5 times

OD measuring, calculation of results

450 nm

Stop Solution

STOP
100 µL

Incubation 2

+18-25 °C
15'

Substrate Solution

SUBS TMB
100 µL

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	3
4. KIT COMPONENTS	4
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	5
6. WARNING AND PRECAUTIONS	5
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	6
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	6
9. REAGENTS PREPARATION	7
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	9
12. EXPECTED VALUES	9
13. PERFORMANCE CHARACTERISTIC	9
14. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
thyroxine in human serum or plasma
T4 EIA

1. INTENDED USE

The T4 EIA kit is an enzyme immunoassay, intended for the quantitative determination of thyroxine in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Thyroxine (T4) and triiodothyronine (T3) are hormones that are produced by the thyroid gland and circulate in the blood both free and bound - mainly with thyroxine-binding globulin (TBG). Only free T3 and T4 are characterized by Hormonal activity, but their share is very small: 0.03% of the total content for T4 and 0.3% - for T3. Concentration of T4 in serum blood is the most accepted indicator of thyroid gland function, which allows you to clearly distinguish between hyper-, hypo- and euthyroidism.

Increase of total T4 concentration is observed with hyperthyroidism, with pituitary tumors, with conditions with elevated TSH levels (pregnancy, acute or chronic active hepatitis, estrogen-secreting tumors or estrogen intake, genetically conditional increase), while taking oral contraceptives, heroin, methadone, thyroid drugs, TSH, thyroliberin.

Decrease of total T4 concentration is observed in hypothyroidism, panhypopituitarism, states of low levels of TSH (acromegaly, nephrotic syndrome, hypoproteinemia, chronic liver disease, androgen-secreting tumors, or androgens, genetically determined decrease), hemolysis, exercise, when taking amino salicylic and acetylsalicylic acids, glucocorticoids, sulfonamides, cholestyramine, reserpine, potassium iodide, triiodothyronine.

3. TEST PRINCIPLE

Determination of the thyroxine is based on competition principle of the enzyme immunoassay. Microwells plate is coated with specific murine monoclonal to thyroxine antibodies. Thyroxine conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage thyroxine from the specimen competes with the conjugated thyroxine for coating antibodies. As a result, a complex bounded to the solid phase and containing peroxidase is formed.
- during the second stage, the complexes formed due the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. Optical density in the microwell is inversely related to the quantity of the measured thyroxine in the specimen of the serum (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of thyroxine in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P212Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to T4; ready to use
C212Z	CAL 1	Calibrator C1	0.5 mL	1	Solution based on human plasma, free of thyroxin, with preservative, ready to use (yellow liquid)
C212Z	CAL 2-5	Calibrators	0.5 mL	4	Solutions based on human plasma, containing 32; 64; 160 and 320 nmol/L of thyroxin, with preservative, ready to use (red liquids)
Q212Z	CONTROL	Control Serum	0.5 mL	1	Solution based on human plasma, containing of known thyroxin content, with preservative, ready to use (colourless liquid)
T212XZ	CONJ	Conjugate Solution	14 mL	1	Solution of thyroxin conjugated to the horseradish peroxidase; ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent; 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)

The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength;
- dry thermostat for $+37^{\circ}\text{C}\pm 2^{\circ}\text{C}$;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 μL ;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The T4 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The T4 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- *NOTE: Single freezing of Calibrators and Control Serum in aliquots is allowed*
- diluted washing solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of washing solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 12 wells for Calibrators and Control Serum (2 wells for each calibrator (CAL 1-5) and 2 wells for control serum (Q)).
- 10.2 Dispense **25 µL of Calibrators and Control Serum as well as 25 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

Note: during performing several independent series of tests, Calibrators, and Control Sample should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP3	SAMP3	SAMP11	SAMP11						
B	CAL2	CAL2	SAMP4	SAMP4	SAMP12	SAMP12						
C	CAL3	CAL3	SAMP5	SAMP5								
D	CAL4	CAL4	SAMP6	SAMP6								
E	CAL5	CAL5	SAMP7	SAMP7								
F	Q	Q	SAMP8	SAMP8								
G	SAMP1	SAMP1	SAMP9	SAMP9								
H	SAMP2	SAMP2	SAMP10	SAMP10								

- 10.3 Add **100 µL of the Conjugate Solution** to all wells.
- 10.4 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C**.
- 10.5 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL.
- 10.6 Add **100 µL of Substrate Solution** to all wells. The introduction of the substrate solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.7 Add **100 µL of Stop Solution** to all wells in the same order as the substrate solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.8 Read the optical density (OD) of the wells at 450nm using a microplate photometer within 5 minutes of adding the Stop Solution.
- 10.9 Plot a calibration curve in semi-logarithmic coordinates: (x) is the decimal logarithm of the T4 concentration in the calibrators nmol/L, (y) – OD versus T4 concentration (OD 450 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve. Adjust the concentration of CAL1 to an infinitesimally small value, for example, 0.001 nmol/L.
- 10.10 Determine the corresponding concentration of T4 in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 1.2, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

12.1. Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for T4. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying

NOTE: values of T4 concentrations in the tested samples that are below the LoD (3.0 nmol/L) and also exceed the value of the upper calibrator (320 nmol/L) should be provided in the following form : «the T4 concentration of tested sample X is «lower than 3.0 nmol/L» or «higher than 320 nmol/L».

12.2. The calibrators concentration values of the T4 EIA kit are expressed in nmol/L. To calculate concentrations in µg/dl, the received concentration value in nmol/L shall be multiplied by 0.0775.

$$1 \text{ nmol/L} = 0.0775 \text{ µg/dl}$$

Sex, age	Units, nmol/L		Units alternative, µg/dl	
	Lower limit	Upper limit	Lower limit	Upper limit
Healthy donors	60	160	4.7	12.4
Males				
>61 yrs	60	129	4.7	10.0
Females				
>61 yrs	70	135	5.4	10.5
Children				
1-5 yrs	90	190	7.0	14.7
6-10 yrs	83	170	6.4	13.2
>10 yrs	60	160	4.7	12.4

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

3.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, nmol/L	CV, %
1	17.5	4.36
2	110.7	3.67

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, nmol/L	CV, %
1	16.4	1.17
2	111.1	5.43

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, nmol/L	Concentration2, nmol/L	Concentration3, nmol/L	CV, %
1	14.59	13.67	15.39	5.92
2	116.23	114.53	120.13	2.45

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known T4 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 0.75–15 nmol/L $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest T4 concentration in the serum or plasma sample that is detected by the T4 EIA kit is no lower than 3 nmol/L.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for T4 EIA kit is 32 nmol/L.

13.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of T4 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
T3	0.5
D-Thyroxin	30

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
free thyroxin in human serum or plasma

fT4 EIA

Catalogue number **REF** **K214**

For 96 determinations

In vitro diagnostic medical device

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Authorized Representative in EU:
Polmed.de Beata Rozwadowska
Fichtenstr. 12A, 90763 Fuerth, Germany
tel.:+ 49 911 931 639 67
E-mail: info@polmed.de
www.polmed.de

ASSAY PROCEDURE

Dispensing of Calibrators, Control Serum and test samples

test samples

SORB MTP

CAL

CAL

CONTROL

25 μ l

Conjugate Solution

CONJ HRP

100 μ l

Incubation 1

+37 $^{\circ}$ C

60'

Washing 5 times

Washing 5 times

OD measuring, calculation of results

450 nm

Stop Solution

STOP

100 μ l

Incubation 2

+18-25 $^{\circ}$ C

15'

Substrate Solution

SUBS TMB

100 μ l

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	2
4. KIT COMPONENTS	3
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	4
6. WARNING AND PRECAUTIONS	4
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	5
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	5
9. REAGENTS PREPARATION	6
10. ASSAY PROCEDURE	6
11. TEST VALIDITY	7
12. EXPECTED VALUES	7
13. PERFORMANCE CHARACTERISTICS	8
14. REFERENCES	9
SAMPLES IDENTIFICATION PLAN	10

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
free thyroxin in human serum or plasma
fT4 EIA

1. INTENDED USE

The fT4 EIA kit is an enzyme immunoassay, intended for the quantitative determination of free thyroxin in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Thyroid hormones thyroxin (T4) and 3,5,3'-triiodothyronine (T3) exert regulatory influences on growth, differentiation, cellular metabolism and development of skeletal and organ systems. T4 and T3 in blood are found both in free and bound form – mostly, they are bound to thyroxin binding globulin (TBG). Only free forms of T3 and T4 exert hormonal activity also their percentage is very low – 0.3% for T3 and 0.03% for T4.

The concentration of T4 is generally accepted as an index of thyroid function which provide enough information to differentiate between hyper-, hypo- and euthyroidism.

Elevation of total T4 is found in hyperthyroidism, in patients with tumours of pituitary gland, in subjects with elevated TBG level (pregnancy, acute or chronic active hepatitis, estrogen-secreting tumours or estrogen intake, hereditary elevation of TBG), in patients taking oral contraceptives, heroin, methadone, thyroid preparations, TSH, thyroliberin.

Low total T4 is found in hypothyroidism, in patients with panhypopituitarism, in subjects with low TBG level (acromegaly, nephritic syndrome, hypoproteinemia, chronic liver diseases, androgen-secreting tumours, hereditary reduction), in patients taking aminosalicilic and acetylsalicilic acids, cholestyramine, reserpine, potassium iodide, triiodothyronine.

3. TEST PRINCIPLE

Determination of free thyroxin is based on competition principle of the enzyme immunoassay. Microwells plate is coated with specific murine monoclonal antibodies to T4. fT4 conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage fT4 from the specimen competes with the conjugated fT4 for coating antibodies. As a result, a complex bounded to the solid phase and containing peroxidase is formed.

- during the second stage, the complexes formed due the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. Optical density in the microwell is inversely related to the quantity of the measured fT4 in the specimen of the serum (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of fT4 in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P214Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to T4; ready to use
C214Z	CAL 1	Calibrator C1	0.5 mL	1	Solution based on human plasma, free of FT4, with preservative, ready to use (yellow liquid)
C214Z	CAL 2-6	Calibrators	0.5 ml	5	Solutions based on human plasma, containing 5; 10; 25, 50 and 100 pmol/L of FT4, with preservative, ready to use (red liquids)
Q214Z	CONTROL	Control Serum	0.5 ml	1	Solution based on human plasma, containing of known FT4 content, with preservative, ready to use (colourless liquid)
T214Z	CONJ HRP	Conjugate Solution	12 ml	1	Solution of FT4 conjugated to the horseradish peroxidase; ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	12 ml	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 ml	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	12 ml	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)

The kit also includes instruction for use, quality control data sheet and plate sealing tape (1 pcs.)

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength;
- dry thermostat for $+37^{\circ}\text{C}\pm 1^{\circ}\text{C}$;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 μL ;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below; do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The ft4 EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The ft4 EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- the remaining strips should be immediately resealed in the bag along with the silica gel, closed with the zip-lock, and stored at +2...+8°C within 2 months
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- *NOTE: Single freezing of Calibrators and Control Serum in aliquots is allowed.*
- diluted washing solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. The remaining strips should be immediately resealed in the bag along with the silica gel and closed with the zip-lock to prevent moisture from affecting the plate's strips.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of washing solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 14 wells for Calibrators and Control Serum (2 wells for each calibrator (CAL 1-6) and 2 wells for control serum (Q)).
- 10.2 Dispense **25 µL of Calibrators and Control Serum as well as 25 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

Note: during performing several independent series of tests, Calibrators, and Control Sample should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP2	SAMP2	SAMP10	SAMP10						
B	CAL2	CAL2	SAMP3	SAMP3	SAMP11	SAMP11						
C	CAL3	CAL3	SAMP4	SAMP4	SAMP12	SAMP12						
D	CAL4	CAL4	SAMP5	SAMP5								
E	CAL5	CAL5	SAMP6	SAMP6								
F	CAL6	CAL6	SAMP7	SAMP7								
G	Q	Q	SAMP8	SAMP8								
H	SAMP1	SAMP1	SAMP9	SAMP9								

- 10.3 Add **100 µL of the Conjugate Solution** to all wells.
- 10.4 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C**.
- 10.5 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 µL of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5µL. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the Washing Solution volume can be increased to 350 µL
- 10.6 Add **100 µL of Substrate Solution** to all wells. The introduction of the substrate solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.7 Add **100 µL of Stop Solution** to all wells in the same order as the substrate solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.8 Read the optical density (OD) of the wells at 450nm using a microplate photometer within 5 minutes of adding the Stop Solution.
- 10.9 Plot a calibration curve in semi-logarithmic coordinates: (x) is the decimal logarithm of the fT4 concentration in the calibrators pmol/L, (y) – OD versus fT4 concentration (OD 450 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve. Adjust the concentration of CAL1 to an infinitesimally small value, for example, 0.001 pmol/L.
- 10.10 Determine the corresponding concentration of fT4 in tested samples from the calibration curve.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 1.2, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for fT4. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE: values of fT4 concentrations in the tested samples that are below the LoD (0.75 pmol/L) and also exceed the value of the upper calibrator (100 pmol/L) should be provided in the following form : «the fT4 concentration of tested sample X is «lower than 0.75 pmol/L» or «higher than 100 pmol/L».

Sex, age	Units, pmol/L	
	Lower limit	Upper limit
Healthy donors		
< 60 yrs	10	25
> 60 yrs	10	21
Pregnancy week		
1st trimester	9	26
2nd trimester	6	21
3rd trimester	6	21

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, pmol/L	CV, %
1	54.4	5.83
2	85.23	3.67

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, pmol/L	CV, %
1	54.36	1.15
2	85.73	3.23

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, pmol/L	Concentration2, pmol/L	Concentration3, pmol/L	CV, %
1	54.59	52.67	60.39	7.19
2	85.23	87.53	85.13	1.58

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined that it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known fT4 concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 5-100 pmol/L $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest fT4 concentration in the serum or plasma sample that is detected by the fT4 EIA kit is no lower than 0.75 pmol/L.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for fT4 EIA kit is 5 pmol/L.

13.1.5 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of fT4 with other analytes is shown in the table:

Analyte	Cross-reactivity, %
L-Thyroxin	100
D-Thyroxin	94
3,3',5'-Triiodo-L-Thyronine (Reverse T3)	86
3,3',5'-Triiodo-L-Thyronine (T3)	3.3
3,3',5'-Triiodo-D-Thyronine	1.8
3,3',5'-Triiodothyropropionic acid	0.6

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 044 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
carcinoembryonic antigen
in human serum or plasma

CEA EIA

Catalogue number **REF** **K224**

For 96 determinations

In vitro diagnostic medical device


XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 050 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Authorized Representative in EU:
Polmed.de Beata Rozwadowska
Fichtenstr. 12A, 90763 Fuerth, Germany
tel.:+ 49 911 931 639 67
E-mail: info@polmed.de
www.polmed.de

ASSAY PROCEDURE

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	3
4. KIT COMPONENTS	4
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	5
6. WARNING AND PRECAUTIONS	5
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	6
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	6
9. REAGENTS PREPARATION	7
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	9
12. EXPECTED VALUES	9
13. PERFORMANCE CHARACTERISTICS	9
14. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
carcinoembryonic antigen
in human serum or plasma
CEA EIA

1. INTENDED USE

The CEA EIA kit is an enzyme immunoassay, intended for the quantitative determination of carcinoembryonic antigen in human serum or plasma.

Quantitative determination of carcinoembryonic antigen (CEA) in blood serum (plasma) is used as an auxiliary method of early diagnosis, monitoring of tumors of the gastrointestinal tract, breast, lung and other adenocarcinomas, as well as assessment of the effectiveness of the therapy for all population groups.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Carcinoembryonic antigen (CEA) represents a family of heavily glycosylated glycoproteins with MW 180–200 kDa which is expressed and secreted by normal human gastrointestinal mucosa. An increase of the level of CEA in the blood serum may precede recurrence of colon or rectal cancer, which is registered in an average of 4-6 months before the development of clinical manifestations of relapse. Although up to 30% of patients with cancer relapse of this localization do not have elevated levels of CEA in the serum, periodic determining the concentration of CEA is important for monitoring patients after surgery - an increase in the level of CEA indicates a recurrence of cancer. Increase the level of CEA is noted in a number of other epithelial tumors, including carcinoma of the breast, stomach, bronchi, pancreas, esophagus, ovaries and endometrium. Determining the content of the CEA is of the greatest importance in the evaluation effectiveness of anticancer therapy (chemo-, radio- or immunotherapy), as well as during follow-up of patients after surgical removal of tumors for the purpose timely detection of relapse. In blood circulation, there are the substances showing high degree of similarity to CEA (NCA, NCA2); this fact requires the use of highly specific anti-CEA reagents. In a present test-system, we use for capturing CEA the monoclonal antibody 3C6. Due to high prevalence of serum CEA elevation in benign diseases (mucosal inflammations), this test system is not recommended for screening for malignant tumours.

3. PRINCIPLE OF THE TEST

The determination of the CEA is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to human CEA. Second antibodies – murine monoclonal antibodies to human CEA conjugated to the horseradish peroxidase is used as enzyme conjugate.

The analysis procedure includes two stages of incubation:

- during the first stage CEA from the specimen is captured by the antibodies coated onto the microwell surface, as well as horseradish peroxidase-conjugated monoclonal antibodies bind to free epitopes of immobilized CEA;

- during the second stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured CEA in the serum specimen (plasma). The concentration is determined according to the calibration graph of the dependence of the optical density on the content of CEA in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P224Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to human CEA; ready to use
C224Z	CAL 1	Calibrator C1	6 mL	1	Solution based on tris buffer (pH 7.2-7.4), free of human CEA, with preservative, ready to use (yellow liquid)
C224Z	CAL 2-6	Calibrators	0.8 mL	5	Solution based on tris buffer (pH 7.2-7.4), containing 2; 4; 8; 32 and 64 ng/mL of CEA, with preservative, ready to use (red liquids)
Q224Z	CONTROL	Control Serum	0.8 mL	1	Solution based on human serum, containing of known CEA content, with preservative, ready to use (colourless liquid)
T224Z	CONJ HRP	Conjugate Solution	14 mL	1	Solution of murine monoclonal antibodies to human CEA conjugated to the horseradish peroxidase; ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)

The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength or 450\620-680 nm;
- thermostat shaker maintaining a speed of 300 rpm and temperature of +37°C ±3°C;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 µL;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLE

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The CEA EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The CEA EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;

NOTE: Single freezing of Calibrators and Control Serum in aliquots is allowed.

- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of Washing Solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

9.4. Samples preparation

If suggested analyte concentration in the sample exceeds the 64 ng/mL, additionally dilute this sample accordingly, using (Calibrator C1). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.

NOTE: in order to obtain reliable results, we recommend to use several successive dilutions of the blood serum (plasma) sample

Do not dilute Control Serum and Calibrators!

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 14 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-6) and 2 wells for Control Serum (Q)).
- 10.2 If necessary, dilute the test samples as described in 9.4.
- 10.3 Dispense **100 µL of Conjugate Solution** to all wells.
- 10.4 Dispense **50 µL of Calibrators and Control Serum as well as 50 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators, and Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP2	SAMP2	SAMP10	SAMP10						
B	CAL2	CAL2	SAMP3	SAMP3	SAMP11	SAMP11						
C	CAL3	CAL3	SAMP4	SAMP4	SAMP12	SAMP12						
D	CAL4	CAL4	SAMP5	SAMP5								
E	CAL5	CAL5	SAMP6	SAMP6								
F	CAL6	CAL6	SAMP7	SAMP7								
G	Q	Q	SAMP8	SAMP8								
H	SAMP1	SAMP1	SAMP9	SAMP9								

- 10.5 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C with continuous shaking 300 rpm**.
- 10.6 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 μ L of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5 μ L. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the wash solution volume can be increased to 350 μ L.
- 10.7 Add **100 μ L of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.8 Add **100 μ L of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.9 Read the optical density (OD) of the wells at 450nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the Stop Solution.
- 10.10 Plot a calibration curve in linear coordinates: (x) is the CEA concentration in the calibrators ng/mL, (y) – OD versus CEA concentration (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.
- 10.11 Determine the corresponding concentration of CEA in tested samples from the calibration curve. In the case of preliminary dilution of the test sample (see 9.4), the obtained result should be multiplied by the dilution factor.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for CEA. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE: values of CEA concentrations in the tested samples that are below the LoD (0.5 ng/mL) and also exceed the value of the upper Calibrator (64 ng/mL) should be provided in the following form : «the CEA concentration of tested sample X is «lower than 0.5 ng/mL» or «higher than 64ng/mL».

Sex, age	Units, ng/mL	
	Lower limit	Upper limit
Non-smokers	-	5
Smokers	-	10

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, ng/mL	CV, %
1	36.64	4.62
2	12.43	7.29

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, ng/mL	CV, %
1	11.56	8.91
2	5.44	10.28

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, ng/mL	Concentration2, ng/mL	Concentration3, ng/mL	CV, %
1	3.67	3.41	3.59	11.6
2	7.51	7.37	7.69	8.7

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined whether it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known CEA concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 2-8 ng/mL $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest CEA concentration in the serum or plasma sample that is detected by the CEA EIA kit is no lower than 0.5 ng/mL.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for CEA EIA kit is 1.0 ng/mL.

13.1.5 Hook Effect

Hook effect is absent for all samples up to reasonably foreseen concentrations 500 ng/mL.

13.1.6 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of CEA with other analytes is shown in the table:

Analyte	Cross-reactivity, %
CA 125	< 0.1
CA 19-9	< 0.1
CA 15-3	< 0.1

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LOT _____ DATE _____

	Manufacturer
	<i>In vitro</i> diagnostic medical device
	Catalogue number
 YYYY-MM	Use-by date
	Batch code
	Temperature limit
	Contains sufficient for <n> tests
	Caution
	Consult instructions for use
	Conformity Marking with technical regulations in Ukraine
	Authorized representative in the European Community/European Union
	CE Conformity Marking

**For any issues related to operation of the kit and technical support,
please contact by telefon number**

+38 044 294-69-78

or write to:

qa@xema.com.ua

XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 050 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
alpha-fetoprotein in human serum or plasma

AFP EIA

Catalogue number **REF** **K225**

For 96 determinations

In vitro diagnostic medical device


XEMA LLC
Akademika Yefremova St. 23
03179, Kyiv, Ukraine
tel.:+38 050 422-62-16
tel.:+38 044 294-69-78
E-mail: qa@xema.com.ua
www.xema.in.ua

Authorized Representative in EU:
Polmed.de Beata Rozwadowska
Fichtenstr. 12A, 90763 Fuerth, Germany
tel.:+ 49 911 931 639 67
E-mail: info@polmed.de
www.polmed.de


ASSAY PROCEDURE

Conjugate Solution


100 μ L

Dispensing of Calibrators, Control Serum and test samples

50 μ L


Incubation 1


+37 $^{\circ}$ C

60'

Washing 5 times


OD measuring, calculation of results


450 /620-680 nm

Stop Solution

100 μ L


Incubation 2

+18-25 $^{\circ}$ C

15'

Substrate Solution

100 μ L

CONTENT

1. INTENDED USE	2
2. GENERAL INFORMATION	2
3. TEST PRINCIPLE	3
4. KIT COMPONENTS	4
5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED	5
6. WARNING AND PRECAUTIONS	5
7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLES	6
8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL	6
9. REAGENTS PREPARATION	7
10. ASSAY PROCEDURE	7
11. TEST VALIDITY	9
12. EXPECTED VALUES	9
13. PERFORMANCE CHARACTERISTICS	9
14. REFERENCES	10
SAMPLES IDENTIFICATION PLAN	11

Instruction for use
A solid-phase enzyme immunoassay kit
for the quantitative determination of
alpha-fetoprotein in human serum or plasma
AFP EIA

1. INTENDED USE

The AFP EIA kit is an enzyme immunoassay, intended for the quantitative determination of alpha-fetoprotein in human serum or plasma.

The field of application is clinical laboratory diagnostics.

2. GENERAL INFORMATION

Alpha-fetoprotein (AFP) is a glycoprotein with a MW ca. 65 kDa which is secreted by fetal liver and yolk sac. AFP represents the main protein of fetal serum while being found in trace quantities in adults. Serum AFP quantitative determination is used in primary diagnostics and monitoring of hepatocellular liver cancer, trophoblastic tumours of testicles and ovary as well as teratomas and teratocarcinomas.

Quantitative determination of AFP in serum of pregnant women or in amniotic fluid during week 15-20 of gestation is widely used for laboratory screening of Down syndrome and defects of spinal cord.

3. PRINCIPLE OF THE TEST

The determination of the alpha-fetoprotein (AFP) is based on the two-site sandwich enzyme immunoassay principle. On the inner surface of the microplate wells are immobilized specific murine monoclonal antibodies to human AFP. Second antibodies – murine monoclonal antibodies to human AFP conjugated to the horseradish peroxidase is used as enzyme conjugate. The analysis procedure includes two stages of incubation:

- during the first stage AFP from the specimen is captured by the antibodies coated onto the microwell surface, as well as horseradish peroxidase-conjugated monoclonal antibodies bind to free epitopes of immobilized AFP;

- during the second stage, the complexes formed due to the reaction with the chromogen 3,3',5,5'-tetramethylbenzidine are visualized.

After stopping the reaction with a stop solution, the intensity of the color of the microwells is measured. The optical density in the microwell is directly related to the quantity of the measured AFP in the serum specimen (plasma).

The concentration is determined according to the calibration graph of the dependence of the optical density on the content of AFP in the calibration samples.

4. KIT COMPONENTS

Code of component	Symbol	Name	Volume	Qty, pcs.	Description
P225Z	SORB MTP	Microplate	-	1	96-well polystyrene strip microplate coated with murine monoclonal antibodies to human AFP; ready to use
C225Z	CAL 1	Calibrator C1	6 mL	1	Solution based on tris buffer (pH 7.2-7.4), free of human AFP, with preservative, ready to use (yellow liquid)
C225Z	CAL 2-6	Calibrators	0.8 mL	5	Solution based on tris buffer (pH 7.2-7.4), containing 5; 15; 50; 150 and 500 IU/mL of AFP, with preservative, ready to use (red liquids)
Q225Z	CONTROL	Control Serum	0.8 mL	1	Solution based on human serum, containing of known AFP content, with preservative, ready to use (colourless liquid)
T225Z	CONJ HRP	Conjugate Solution	14 mL	1	Solution of murine monoclonal antibodies to human AFP conjugated to the horseradish peroxidase, ready to use (red liquid)
R055Z	SUBS TMB	Substrate Solution	14 mL	1	Tetramethylbenzidine (TMB) substrate solution; ready to use (colourless liquid)
S008Z	BUF WASH 26X	26x Concentrate Washing Solution	22 mL	1	Buffer solution with detergent, 26x concentrate (colourless liquid)
R050Z	STOP	Stop Solution	14 mL	1	5.0% solution of sulphuric acid; ready to use (colourless liquid)
The kit also includes instruction for use, quality control data sheet and plate sealing tape (2 pcs.)					

5. EQUIPMENT AND MATERIAL REQUIRED BUT NOT PROVIDED

- microplate photometer with 450 nm wavelength;
- dry thermostat for $+37^{\circ}\text{C}\pm 2^{\circ}\text{C}$;
- automatic plate washer (optional);
- micropipettes with variable volume, range volume 5-1000 μL ;
- graduated cylinder of 1000 mL capacity;
- distilled or deionized water;
- timer;
- vortex mixer;
- disposable gloves;
- absorbent paper.

6. WARNING AND PRECAUTIONS

In order to prevent incorrect results, strictly follow the recommended order and duration of the analysis procedure.

6.1. The kit is for *in vitro* diagnostic use only. For professional laboratory use.

6.2. Follow the rules mentioned below during the kit using:

- do not use kit beyond expire date;
- do not use the kit if its packaging is damaged;
- in order to avoid contamination, use new tips to pipette samples and reagents;
- use only verified equipment;
- close each vial with its own cap, after using the reagent;
- do not use components of other kits or reagents of other manufacturers;
- do not let wells dry after completing the rinsing step; immediately proceed to the next stage;
- avoid bubbles when adding reagents.

ATTENTION! The TMB substrate solution is light sensitive. Avoid prolonged exposure of the component to light.

6.3. Some kit components, such as stop solution, substrate solution, and washing solution, may cause toxic or irritant effects. If they get on the skin or mucosa, the affected area should be washed with plenty of running water.

6.4. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

6.5. The Calibrators and Control Serum included in the kit are negative for antibodies to HIV 1,2, hepatitis C virus and HBsAg, but the reagents should be considered as potentially infectious material and handled carefully.

6.6. Specimens must not contain any azide compounds, as they inhibit activity of peroxidase.

6.7. Wear protective gloves, protective clothing, eye protection, face protection.

6.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

6.9. Safety Data Sheet for this product is available upon request directly from XEMA LLC.

6.10. Serious incidents related to the kit must be reported to the manufacturer, Authorized Representative, and to the Competent Authority of the EU member state(s) where the incident has occurred.

7. SPECIMEN COLLECTION, TRANSPORTATION AND STORAGE OF SAMPLE

7.1. Blood sampling should be carried out from the cubital vein with a disposable needle using a vacuum blood sampling system. Serum or plasma specimens should be clearly labeled and identified. Serum must be separated from the clot as early as possible to avoid hemolysis of red blood cells. If there are any visible particles in the sample, they should be removed by centrifugation at 3000-5000 rpm for 20 minutes at room temperature or by filtration.

Don't use samples with high lipidemia, hemolysis as they may give false test results.

7.2. Specimen should be stored at +2...+8°C up to 3 days. Specimen held for a longer time, should be placed in a freezer at -15°C or below, do not refreeze/thaw samples.

7.3. For the transportation of samples, it is recommended to use triple packaging. The primary package is the labeled tube containing the sample. Secondary packaging is a polyethylene bag that is hermetically closed with a zip-lock. The outer packaging is a heat-insulating container, while the secondary packaging is placed in the outer packaging for transportation in the center of the thermal container. Frozen refrigerants are placed on the bottom, along the side walls of the thermal container, and cover the samples with them.

8. TRANSPORTATION AND STORAGE TERMS OF KIT, WASTE DISPOSAL

Information about the singularity storage conditions, transportation of the kit, and disposal of waste should be taken into account by all persons who participate in these processes.

8.1. Transportation

The AFP EIA kit should be transported in the manufacturer's packaging at +2...+8°C. Single transportation at the temperature up to 25°C for 5 days is acceptable.

8.2. Storage

The AFP EIA kit should be stored in the manufacturer's packaging at +2...+8°C. Do not freeze.

The kit contains reagents sufficient for 96 determinations including Calibrators and Control Serum.

Once opened test-kit is stable for 2 months when stored properly as intended by manufacturer at 2-8°C.

In case of partial use of the kit, the components should be stored in the following way:

- strips that remain unused must be carefully sealed with the plate sealing tape and stored at +2...+8°C within 2 months;
- Substrate Solution, Stop Solution, and Washing Solution concentrate after opening the vial, can be stored tightly closed at +2...+8°C until the kit's shelf life;
- Conjugate Solution, Calibrators and Control Serum after opening the vial, can be stored tightly closed at +2...+8°C within 2 months;
- *NOTE: Single freezing of Calibrators and Control Serum in aliquots is allowed.*
- diluted Washing Solution can be stored at room temperature (+18...+25°C) for up to 5 days or at +2...+8°C for up to 14 days.

Kits that were stored in violation of the storage condition cannot be used.

8.3. Disposal

Expired kit components, used reagents and materials, as well as residual samples must be inactivated and disposed of in accordance with legal requirements.

9. REAGENTS PREPARATION

9.1. All reagents (including microstrips) and test samples should be allowed to reach room temperature (+18...+25 °C) for at least 30 minutes before use.

9.2. Microplate preparation

Open the package with the microplate and install the required number of strips into the frame. Unused strips must be sealed with plate sealing tape to prevent moisture from affecting the plate's holes and placed back in the bag.

9.3. Washing Solution preparation

Add the contents of the 22 mL Washing Solution concentrate vial to 550 mL of distilled or deionized water and mix thoroughly. In case of partial use of the kit, take the necessary amount of Washing Solution concentrate and dilute it 26 times with distilled or deionized water.

The spending of the components in case of partial use of the kit is given in the table:

Quantity of strips	1	2	3	4	5	6	7	8	9	10	11	12
Volume of the Washing Solution concentrate, mL	1.8	3.6	5.4	7.2	9	10.8	12.6	14.4	16.2	18	19.8	22
Volume of water, mL	45	90	135	180	225	270	315	360	405	450	495	550

9.4. Samples preparation

If suggested analyte concentration in the sample exceeds the 500 IU/mL, additionally dilute this sample accordingly, using (Calibrator C1). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.

NOTE: in order to obtain reliable results, we recommend to use several successive dilutions of the blood serum (plasma) sample

Do not dilute Control Serum and Calibrators!

10. ASSAY PROCEDURE

- 10.1 Put the desired number of strips into the frame based on the number of test samples in 2 replicates and 14 wells for Calibrators and Control Serum (2 wells for each Calibrator (CAL 1-6) and 2 wells for Control Serum (Q)).
- 10.2 If necessary, dilute the test samples as described in 9.4.
- 10.3 Dispense **100 µL of Conjugate Solution** to all wells.
- 10.4 Dispense **50 µL of Calibrators and Control Serum as well as 50 µL of test serum/plasma samples (SAMP)** to the wells of the microplate according to the scheme below. The introduction of Calibrators, Control Serum and test samples should be carried out within 5 minutes to ensure equal incubation time for the first and last samples.

NOTE: during performing several independent series of tests, Calibrators, and Control Serum should be used each time.

Scheme of introduction of samples

	1	2	3	4	5	6	7	8	9	10	11	12
A	CAL1	CAL1	SAMP2	SAMP2	SAMP10	SAMP10						
B	CAL2	CAL2	SAMP3	SAMP3	SAMP11	SAMP11						
C	CAL3	CAL3	SAMP4	SAMP4	SAMP12	SAMP12						
D	CAL4	CAL4	SAMP5	SAMP5								
E	CAL5	CAL5	SAMP6	SAMP6								
F	CAL6	CAL6	SAMP7	SAMP7								
G	Q	Q	SAMP8	SAMP8								
H	SAMP1	SAMP1	SAMP9	SAMP9								

- 10.5 Carefully mix the contents of the microplate in a circular motion on a horizontal surface, cover strips with a plate sealing tape and incubate for **60 minutes at +37°C**.
- 10.6 At the end of the incubation period, remove and discard the plate cover. Aspirate and wash each well 5 times using an automatic washer or an 8-channel dispenser. For each washing, add 300 μ L of Washing Solution (see 9.3) to all wells, then remove the liquid by aspiration or decantation. The residual volume of the Washing Solution after each aspiration or decantation should be no more than 5 μ L. After washing, carefully remove the remaining liquid from the wells on the absorbent paper. For the automatic washer/analyzer, the wash solution volume can be increased to 350 μ L.
- 10.7 Add **100 μ L of Substrate Solution** to all wells. The introduction of the Substrate Solution into the wells must be carried out within 2-3 minutes. Incubate the microplate in the dark **at room temperature (+18...+25°C) for 15 minutes**.
- 10.8 Add **100 μ L of Stop Solution** to all wells in the same order as the Substrate Solution. After adding the Stop Solution, the contents of the wells turn yellow.
- 10.9 Read the optical density (OD) of the wells at 450nm and reference light filters 620–680 nm using a microplate photometer within 5 minutes of adding the Stop Solution. Set photometer blank on CAL1.
- 10.10 Plot a calibration curve in linear coordinates: (x) is the AFP concentration in the calibrators IU/mL, (y) – OD versus AFP concentration (OD 450 nm / 620–680 nm). Manual or computerized data reduction is applicable at this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve
- 10.11 Determine the corresponding concentration of AFP in tested samples from the calibration curve. In the case of preliminary dilution of the test sample (see 9.4), the obtained result should be multiplied by the dilution factor.

11. TEST VALIDITY

The test run shall be considered valid if the OD of CAL1 is above 0.15, and the values of the Control Serum fall into the required range (see Quality control Data Sheet).

12. EXPECTED VALUES

12.1. Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for AFP. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

NOTE: values of AFP concentrations in the tested samples that are below the LoD (0.9 IU/mL) and also exceed the value of the upper Calibrator (500 IU/mL) should be provided in the following form : «the AFP concentration of tested sample X is «lower than 0.9 IU/mL» or «higher than 500IU/mL».

12.2. The calibrators concentration values of the AFP EIA kit are expressed in IU/mL. To calculate concentrations in ng/mL, the received concentration value in IU/mL shall be multiplied by 1.25.

$$1 \text{ IU/mL} = 1,25 \text{ ng/mL}$$

Sex, age	Units, IU/mL		Units alternative, ng/mL	
	Lower limit	Upper limit	Lower limit	Upper limit
Healthy donors	-	10.0	-	12.5

Medians and SKO (recommended normal range 0.5-2.0)

Pregnancy, week	Median, IU/mL	SKO
14	21.7	0.43
15	28.3	0.47
16	30.0	0.51
17	36.3	0.49
18	43.8	0.52
19	53.3	0.50
20	60.0	0.55
21	63.3	0.57

13. PERFORMANCE CHARACTERISTICS

13.1. Analytical performance characteristics

13.1.1 Precision of Measurement

Repeatability (Intra assay repeatability) was determined by evaluation the coefficient of variation (CV) for 2 different samples during 1 day in 24 replicates on one series of ELISA kit.

Sample	Concentration, IU/mL	CV, %
1	54.16	1.83
2	289	7.67

Reproducibility (Inter assay reproducibility) was determined by evaluating the coefficients of variation for 2 samples during 5 days in 8-replicate determinations.

Sample	Concentration, IU/mL	CV, %
1	55.34	5.75
2	288.47	2.63

Reproducibility between lots was investigated by testing samples for one day on three lots. Each sample was run in 8 replicates.

Sample	Concentration1, IU/mL	Concentration2, IU/mL	Concentration3, IU/mL	CV, %
1	54.6	55.89	54.89	2.95
2	289.4	281.75	283.46	1.41

13.1.2 Trueness

The trueness of measurement is the degree of closeness of the average value obtained from a large number of measurement results to the true value. The bias of the measurement result (bias of measurements) is the difference between the mathematical expectation of the measurement result and the true value of the measurand. The bias was calculated for each sample and it was determined whether it corresponds to the specified limits of $\pm 10\%$.

13.1.3 Linearity

Linearity was determined using sera samples with known AFP concentration (low and high) and mixing them with each other and buffer solution in different proportions. According to the measurements, linear range of kit is 5-500 IU/mL $\pm 10\%$.

13.1.4 Analytical sensitivity

Limit of detection (LoD) – the lowest AFP concentration in the serum or plasma sample that is detected by the AFP EIA kit is no lower than 0.9 IU/mL.

Limit of quantification (LoQ) – the lowest concentration of the analyte in the sample that is determined quantitatively with the declared trueness for AFP EIA kit is 5.0 IU/mL.

13.1.5 Hook Effect

Hook effect is absent for all samples up to reasonably foreseen concentrations 12000 IU/mL.

13.1.6 Analytical specificity

For the analysis result is not affected by the presence in the sample of bilirubin in a concentration of up to 0.21 mg/mL and hemoglobin in a concentration of up to 10 mg/mL.

The cross-reactivity of AFP with other analytes is shown in the table:

Analyte	Cross-reactivity, %
Albumin	< 0.1
hCG	< 0.1
Lactogen	< 0.1

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