## BILIRUBIN TOTAL & DIRECT 4+1

### Referințe:

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### BIBLIOGRAFIE

- Higgins, T., et al., *Hemoglobin, Iron, and Bilirubin*, Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds, Philadelphia USA), (2008), 509.
- Dufour, D.R. *The liver: Function and chemical Pathology*, Clinical Chemistry: Theory, Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc.), (2010), 586 and appendix.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
- Wu, H.B., *General Clinical Tests*. Tietz Clinical guide to laboratory tests, 4<sup>th</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 172.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).8.
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
- Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> edition, AACC Press (1997).
- Young D.S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> edition, AACC Press (1995).

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## English - EN

### INTENDED USE

This *in vitro* diagnostic device is intended for the calibration of ELITech Clinical Systems quantitative tests listed in the value sheet.  
This *in vitro* diagnostic device is for professional use only.

### COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.  
- Sodium azide < 0.1 % (w/w)  
- Concentrations for each analyte to test are lot-specific.

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### TRACEABILITY

The traceability is indicated in the value sheet enclosed in the kit.

### PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved or cleared in compliance with European Directive 98/79/EC, Annex II, List A. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.  
- Take precautions when handling broken glass vials as sharp edges can injure the user.  
- These products contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.  
- Take normal precautions and adhere to good laboratory practice.  
- Use clean or single use laboratory equipment only to avoid contaminations.  
- Consult Safety Data Sheet (SDS) for a proper handling.

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### PREPARATION

- Carefully open the vial, avoiding the loss of lyophilizate.  
- Add in exactly 3 mL of distilled or deionized water.  
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

### PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particules would indicate deterioration.  
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).  
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

### STABILITY

**Prior to reconstitution :**  
- Stored at 2-8 °C and protected from light.  
- Do not use after expiration date indicated on the vial label.

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- Stability of the components:  
Between 15 and 25 °C : 8 hours  
Between 2 and 8 °C : 2 days  
Between -25 and -15 °C : 4 weeks (when frozen once)

### Exceptions:

- Stability of total bilirubin (when stored protected from light):  
Between 15 and 25 °C : 6 hours  
Between 2 and 8 °C : 1 day  
Between -25 and -15 °C : 2 weeks (When frozen once)  
- Stability of direct bilirubin (when stored protected from light):  
Between 15 and 25 °C : 3 hours  
Between 2 and 8 °C : 8 hours  
Between -25 and -15 °C : 2 weeks (when frozen once)

**Note:** Store vials tightly capped and protected from light.

### PROCEDURE

To use ELICAL 2, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

### LIMITATIONS

Using ELICAL 2 has been validated with the ELITech systems (Analyzers and reagents listed in the value sheet).

Using any other system should be validated by the laboratory.

### VALUES

The concentrations are indicated in the value sheet enclosed in the kit.  
For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

### DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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## English - EN

### INTENDED USE

These *in vitro* diagnostic devices are intended for the quality control of the performances of ELITech Clinical Systems quantitative tests listed in the value sheet. These *in vitro* diagnostic devices are for professional use only.

### COMPOSITION

- Lyophilized product prepared from human serum spiked with chemical and biological additives.  
- Concentrations for each analyte to test are lot-specific.

### MATERIALS REQUIRED BUT NOT PROVIDED

- ELITech Clinical Systems reagents listed in the value sheet.  
- General Laboratory equipment.

### TRACEABILITY

The traceability is indicated on the instructions for use of the ELITech Clinical Systems reagents used in combination with the recommended calibrator.

### PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.  
- Take precautions when handling broken glass vials as sharp edges can injure the user.  
- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contaminations.

- Consult Safety Data Sheet (SDS) for a proper handling.

### WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

### PREPARATION

- Gently tap the vial so that the control material can be at the bottom of the vial.  
- Carefully open the vial, avoiding the loss of lyophilizate.  
- Add exactly 5 mL of distilled or deionized water.  
- Carefully close the vial and dissolve the contents completely within 30 minutes by occasional gentle stirring avoiding the formation of foam.

**Important:** To reactivate the alkaline phosphatase, allow the reconstituted control serum to stand for one hour at 15-25 °C.

### PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particules would indicate deterioration.  
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).

- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

### STABILITY

#### Prior to reconstitution :

- Store at 2-8 °C and protect from light. Do not freeze.

- Do not use after the expiry date stated on the vial label.

#### After reconstitution :

- These products should be immediately and tightly capped to prevent contamination and evaporation.

- Stability of the components:

Between 2 and 8 °C : 5 days

Between -25 and -15 °C : 4 weeks (when frozen once)

#### Exceptions:

- Stability of total bilirubin (when stored protected from light) :  
Between 2 and 8 °C : 24 hours  
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of direct bilirubin (when stored protected from light) :  
Between 2 and 8 °C : 8 hours  
Between -25 and -15 °C : 5 days (when frozen once)

- Stability of alkaline phosphatase:

Between 2 and 8 °C : 24 hours

Between -25 and -15 °C : 4 weeks (when frozen once)

#### Note:

- Store the vials tightly capped and protected from light after reconstitution.

- Mix content thoroughly prior to use.

- Thawing frozen specimen at 37°C, and vigorous mixing (vortex), should be avoided as this may denature ferritin.

### PROCEDURE

To use ELITROL I & II, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

### LIMITATIONS

Using ELITROL I & II has been validated with the ELITech systems (Analyzers and reagents listed in the value sheet).  
Using any other system should be validated by the laboratory.

### VALUES

The concentrations are indicated in the value sheet enclosed in the kit.  
Individual laboratories may obtain values different from those announced. Technique, equipment and experimental error may produce slightly different values. Each laboratory should determine their own mean values.

For users of Selectra instruments allowing data import for tests, calibrators and controls, use corresponding barcode, available on the value sheet.

### DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

### ASSISTANCE TECHNIQUE

Contact your local distributor or ELITech Clinical Systems SAS (CCsupport@elitechgroup.com).

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
☞ **BIBLIOGRAFIE**

1. Wu, A.H.B., Clinical guide to laboratory tests, 4<sup>th</sup> Ed., (W.B. Saunders eds.), (2006), 66.
2. Dufour D. R., The Liver: Function and Chemical Pathology. Clinical Chemistry: Theory, Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 586 and appendix.
3. Johnson, A. M., Amino Acids and Proteins. Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 286.
4. Doumas, B. T., et al., Albumin standards and the measurement of serum albumin with bromocresol green, Clin Chim Acta (1971), **31**, 87.
5. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
6. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
7. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
8. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
9. Interference Testing in Clinical Chemistry ; Approved Guideline—Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25**(27).
10. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
11. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> edition, AACC Press (1997).
12. Young D.S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> edition, AACC Press (1995).

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### BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzyme, Tietz Fundamentals of Clinical Chemistry*, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4<sup>th</sup> Ed., (W.B. Saunders Company), (2006), 154
3. Dufour, R., *The liver: Function and chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 586, appendix
4. Ward, M.K., Cockayne, S., *Enzymology. Clinical Chemistry : Concepts and Application*, Anderson, S.C., Cockayne, S. (W.B. Saunders eds. Philadelphia USA), (1993), 238.
5. Schumann, G., et al *Clin Chem Lab Med.*, (2002), **40**, 725.
6. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations*, WHO/DIL/LAB/99.1 Rev.2, (2002), 21.
7. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
8. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
9. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
10. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

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12. Vassault A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686

13. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry ana-lysis of blood samples containing monoclonal immuno-globulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.

14. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).

15. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).

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### ☞ BIBLIOGRAFIE

1. Wu, A.H.B., Tietz Clinical guide to laboratory test, 4<sup>th</sup> Ed., (W.B. Saunders Company), (2006), 684.
2. Itani, O., Tsang, R.C., Bone disease, Clinical Chemistry: Theory Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 614 and appendix.
3. Foley, K. F., Boccuzzi, L., Urine Calcium: Laboratory Measurement and Clinical Utility, Labmedicine, (2010), **41**, 683.
4. Bauer, P. J., Anal. Biochem., (1981), 110, 61.
5. Guder, W.G., et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Endres, D.B., Rude, R. K., Disorders of Bone, Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 711.
7. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
8. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
9. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
10. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
11. Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25**(27).
12. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
13. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> edition, AACC Press (1997).
14. Young D.S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> edition, AACC Press (1995).

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



FTRO-CRCO-v23(12/2018)\_PIT-CRCO-4-v23

## BIBLIOGRAFIE

1. Allston, C.A., *Non protein nitrogenous compounds and renal function*. Clinical Chemistry: Concepts and Application, Anderson, S.C., Cockayne, S. (W.B. Saunders eds. Philadelphia USA), (1993), 369.
2. Newman, D.J., Price C.P., *Non protein nitrogen metabolite*. Tietz Fundamentals of Clinical Chemistry, 5<sup>th</sup> Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 414.
3. Butler, A.R., The Jaffe reaction. *Identification of the coloured species*. Clin. Chim. Acta., (1975), 59, 227.
4. Vasiliades, J., *Reaction of alkaline picrate with creatinine. 1. Kinetics and mechanism of formation of the mono-creatinine picric acid complex*. Clin. Chem., (1976), **22**, 1664.
5. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
6. Tietz, N.W. *Clinical guide to laboratory tests*, 3<sup>rd</sup> Ed, (W.B. Saunders eds. Philadelphia USA), (1995), 186.
7. Vassault A., *et al.*, Ann. Biol. Clin., (1986), 44, 686.
8. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
9. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> edition, AACC Press (1997).
10. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> edition, AACC Press (1995).

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
FTRO-PASL-v21(12/2018)\_PIT-PASL-4-v21

**BIBLIOGRAFIE**

- Henderson, A.R., Donald W.M., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 5<sup>th</sup> Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 352.
- Scherwin, J.E., *Liver fonction. Clinical Chemistry: Theory, Analysis, Correlation*, 4<sup>th</sup> Ed., Kaplan, L.A, Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 492 and appendix.
- German Society for Clinical Chemistry, *Z. Klin. Chem. Klin. Biochem.*, (1972), **10**, 281.
- Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology, *Scand. J. Clin. Lab. Invest.*, (1974), **33**, 291.
- Guder, World Health Organization, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*, WHO/ DIL/LAB/99.1 Rev.2, (2002), 21.
- Rosalki, S.B., *Clin. Chem.*, (1993), **37**, 648.
- Vassault A., *et al.*, *Ann. Biol. Clin.*, (1986), **44**, 686.
- Vassault A., *et al.*, *Ann. Biol. Clin.*, (1999), **57**, 685.
- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immuno-globulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).

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#### BIBLIOGRAFIE

1. Naito, H.K., *Coronary Artery Disease and Disorders of Lipid Metabolism*. *Clinical Chemistry: Theory, Analysis, Correlation*, 4<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C. (Mosby, Inc. eds. St Louis USA), (2003), 603.
2. Expert Panel on Detection, *Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III)*, *Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP)*. *JAMA*, (2001), **285**, 2486.
3. Fossati, P., Prencipe, L., *Serum triglycerides determined colorimetrically with an enzyme that produces hydrogen peroxide*. *Clin. Chem.*, (1982), **28**, 2077.
4. Tietz, N.W., *Clinical guide to laboratory tests*, 3<sup>rd</sup> Ed. (W.B. Saunders eds. Philadelphia USA), (1995), 610.
5. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), **23** (16).

6. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).

7. *Method Comparison and Bias estimation Using Patient Samples; Approved Guideline-Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

8. *Interference Testing in Clinical Chemistry ; Approved Guideline-Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).

9. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.

10. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> edition, AACC Press (1997).

11. Young D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> edition, AACC Press (1995).

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**BIBLIOGRAFIE**

1. Wu, A.H.B., Tietz Clinical guide to laboratory test, 4<sup>th</sup> Ed., (W.B. Saunders Company), (2006), 470
2. Panteghini, M., Bais, R., Enzyme, Tietz Fundamentals of Clinical Chemistry, 6<sup>th</sup> Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
3. Dufour, R., The liver function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation, 5<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 586, appendix.
4. Schumann, G., et al., IFCC Primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 6. Reference procedure for the measurement of catalytic concentration of  $\gamma$ -Glutamyltransferase. Clin. Chem. Lab. Med., (2002), **40**, 734.
5. Guder, W.G., et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples, WHO/DIL/LAB/99.1 Rev.2, (2002).
6. Kytzia H-J, Reference intervals for GGT according to the new IFCC 37°C reference procedure, Clin. Chem. Lab. Med., (2005), **43**, A69.
7. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
8. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
9. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
10. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

11. Interference Testing in Clinical Chemistry ; Approved Guideline—Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).

12. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.

13. Young, D. S., Effects of preanalytical variables on clinical laboratory tests, 2<sup>nd</sup> Ed., AACC Press, (1997).

14. Young, D. S., Effects of drugs on clinical laboratory tests, 4<sup>th</sup> Ed., AACC Press, (1995).

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**BIBLIOGRAFIE**

1. Scherwin, J.E. *Liver function. Clinical Chemistry: Theory, Analysis, Correlation*, 4<sup>th</sup> Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 492 and appendix.
2. Tietz, N.W., *Clinical guide to laboratory tests*, 3<sup>rd</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 518.
3. Doumas, B.T., et al., *Clin. Chem.*, (1981), **27**, 1642.
4. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
5. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI (NCCLS) document EP6-A (2003), **23** (16).
6. *Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
7. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
8. *Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
9. *Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
10. Vassault, A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686.
11. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
12. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACC Press, (1997).
13. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACC Press, (1995).

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