Medtronic

IN.PACT™ Admiral™

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
English

1 Description

The IN.PACT Admiral drug-coated balloon catheter (DCB) is a paclitaxel-coated, over-the-wire (OTW) peripheral balloon catheter manufactured by Medtronic, Inc. The FreePac™ drug coating on the balloon of the IN.PACT Admiral DCB consists of the drug paclitaxel and the excipient urea. Urea facilitates the release and transfer of paclitaxel into the arterial vessel wall. The balloon catheter physically dilates the vessel lumen by percutaneous transluminal angioplasty, and the drug coating is intended to reduce the proliferative response that is associated with restenosis. Paclitaxel stabilizes microtubules to reduce cell proliferation.

The IN.PACT Admiral DCB has a dual-lumen shaft. This shaft branches into two tubes at the proximal end; one tube forms the entrance to the central lumen for the guidewire, while the other tube serves as the passage for the mixture of contrast medium and saline solution that inflates and deflates the balloon. The catheter construction and the balloon material are designed to reach targeted balloon diameters, which will vary depending on the balloon size and defined pressure. The length of each balloon is specified. For all balloon lengths, 2 radiopaque platinum-iridium markers indicate the working length of the balloon to aid in positioning the balloon across the target lesion during fluoroscopy; the 200 mm and 250 mm balloons have 2 additional centrally-placed markers. Guidewire and introducer sheath compatibility is listed in *Table 1*.

Table 1. Compatibility

Minimum introducer sheath compatibility	
Balloon diameter (mm)	Introducer sheath compatibility (Fr)
4.0	5
5.0	6
6.0	
7.0	7
8.0	
9.0	

Table 1. Compatibility (continued)

Minimum introducer sheath compatibility		
10.0	12.0	9
Guidewire compatibility		

The catheter is compatible with a guidewire diameter of 0.035 in (0.89 mm)

The IN.PACT Admiral DCB is available in an assortment of balloon sizes. Nominal balloon diameter and lengths are printed on the hub. The IN.PACT Admiral DCB is not made with natural rubber latex; however, during the manufacturing process, it may have incidental contact with latex.

2 Device composition

The catheter is composed of polyamide (PA12) and polyether block amide (PEBA), with a polycarbonate hub, polyolefin strain relief and platinum/iridium radiopaque markers.

The balloon is composed of polyamide (PA12).

The drug content in the FreePac coating is defined in *Table 2*.

Table 2. Product matrix and paclitaxel content

Balloon diameter (mm) x balloon length (mm)	Nominal paclitaxel content (µg)	Nominal urea content (µg)
4 x 40	1969	276
4 x 60	2848	399
4 x 80	3728	522
4 x 120	5487	768
4 x 150	6807	953
4 x 200	9006	1261
4 x 250	11205	1569
5 x 40	2553	357
5 x 60	3653	511
5 x 80	4752	665
5 x 120	6951	973
5 x 150	8601	1204
5 x 200	11350	1589
5 x 250	14098	1974
6 x 40	3170	444
6 x 60	4489	628
6 x 80	5809	813
6 x 120	8448	1183
6 x 150	10427	1460
6 x 200	13726	1922
6 x 250	17024	2383
7 x 40	3819	535
7 x 60	5358	750
7 x 80	6897	966
8 x 40	4494	629
8 x 60	6253	875
8 x 80	8012	1122
9 x 40	5204	729
9 x 60	7183	1006
9 x 80	9162	1283
10 x 40	5943	832
12 x 40	7522	1053

3 Intended purpose

The intended purpose of the IN.PACT Admiral DCB is to restore lumen patency and blood flow to peripheral arteries, previously stented peripheral arteries with in-stent restenosis, and native arteriovenous dialysis fistulae and to reduce the proliferative response associated with restenosis.

3.1 Intended patient population

The IN.PACT Admiral DCB is intended for patients in need of percutaneous transluminal angioplasty to improve and maintain luminal diameter in peripheral arteries, previously stented peripheral arteries with in-stent restenosis, and native arteriovenous dialysis fistulae.

This device has not been tested on a pediatric population.

3.2 Indications for use

The IN.PACT Admiral DCB is indicated for percutaneous transluminal angioplasty (PTA) in patients with de novo, restenotic, or in-stent restenotic lesions with lengths up to 360 mm in superficial femoral or popliteal arteries.

The IN.PACT Admiral DCB with balloon lengths ≤ 150 mm is also indicated for patients with obstructive lesions of native arteriovenous dialysis fistulae.

3.3 Contraindications

IN.PACT Admiral DCB is contraindicated for use in:

- Coronary arteries and supra-aortic/cerebrovascular arteries
- Lesions that cannot be crossed with a guidewire
- Pregnant or breast-feeding women
- Patients with known allergies or hypersensitivities to paclitaxel

3.4 Clinical benefits

Note: After the EUDAMED website is launched, the Summary of Safety and Clinical Performance (SSCP) can be found at <https://ec.europa.eu/tools/eudamed> using the Basic UDI-DI: 0763000B00004367.

The IN.PACT Admiral DCB is intended to restore lumen patency and blood flow to peripheral arteries, previously stented peripheral arteries with in-stent restenosis, and native arteriovenous dialysis fistulae and to reduce the proliferative response associated with restenosis.

The IN.PACT Admiral DCB balloon device component physically dilates the vessel lumen by PTA (primary mode of action), and the paclitaxel medicinal substance provides a pharmacological agent targeted towards reducing the injury response that leads to restenosis (secondary mode of action).

The positive impact of the IN.PACT Admiral DCB on patients eligible for PTA in peripheral arteries, previously stented peripheral arteries with in-stent restenosis, and native arteriovenous dialysis fistulae is improved or increased lumen diameter and patency.

Percutaneous transluminal angioplasty (PTA) for obstructive disease of the peripheral arteries (not ISR) – The primary clinical benefits include restoring lumen patency and blood flow and inhibition of restenosis. The measurable patient-relevant outcomes and quantitative clinical data from the Medtronic sponsored IN.PACT SFA trial at 5 years are shown in *Table 3*.

Table 3. IN.PACT SFA trial (5 years)

Parameter	IN.PACT SFA/IN.PACT DCB (N=220 subjects) (N=311 devices)
Primary patency through 1080 days ^a	69.5%
Primary sustained clinical improvement at 36 months	68.7% (114/166)
Freedom from clinically-driven TLR at 1800 days ^a	74.5%
Acute success	
Device success (per device)	99.0% (308/311)
Procedural success (per subject)	99.5% (219/220)
Clinical success (per subject)	99.1% (218/220)
Cumulative complications within 1800 days	
MAE composite (death, major target limb amputation, CD-TVR, thrombosis)	42.9% (79/184)
• Death (all-cause)	15.8% (29/184)
• CD-TVR	29.3% (54/184)
• Major target limb amputation	0.5% (1/184)
• Thrombosis	2.2% (4/184)
CD-TLR	25.5% (47/184)
Any TVR	29.9% (55/184)
Any TLR	26.6% (49/184)
Other major secondary endpoints	
Change in quality of life from baseline by EQ-5D index to 36 months	
N	156
Mean ± SD	0.0832 ± 0.2293
Median	0.0905
Min, Max	-0.759, 0.693
Walking impairment by WIQ (%) to 36 months	
N	158
Mean ± SD	71.8 ± 34.2
Median	100.0
Min, Max	0, 100
Nights in hospital due to index lesion through 60 months	
N	219
Mean ± SD	1.6 ± 2.8
Median	1.0
Min, Max	0, 26
Time to first clinically-driven TLR (days) through 60 months	
N	47
Mean ± SD	807.5 ± 433.9
Median	770.0
Min, Max	1, 1701

^a Numbers are cumulative incidence % (number of failures) based on the Kaplan-Meier method.**Table 4.** IN.PACT global/clinical cohort (5 years)

Parameter	IN.PACT global clinical cohort/IN.PACT DCB (N=1406 subjects) (N=3006 devices)
Freedom from clinically-driven TLR at 1800 days ^a	69.4%
Acute success	
Device success (per device)	99.4% (2988/3006)
Procedural success (per subject)	99.3% (1386/1396)
Clinical success (per subject)	98.6% (1376/1396)
Cumulative complications within 1800 days	
MAE composite (death, major target limb amputation, CD-TVR, thrombosis)	48.5% (589/1215)
• Death (all-cause)	20.1% (244/1215)
• CD-TVR	31.4% (381/1215)
• Major target limb amputation	1.6% (19/1215)
• Thrombosis	6.0% (73/1215)
CD-TLR	30.1% (366/1215)
Any TVR	32.1% (390/1215)
Any TLR	30.8% (374/1215)
Other major secondary endpoints	
Change in quality of life from baseline by EQ-5D index to 36 months	
N	823
Mean ± SD	0.1320 ± 0.3328
Median	0.1050
Min, Max	-0.912, 1.150
Walking impairment by WIQ (%) to 36 months	
N	831
Mean ± SD	74.4 ± 30.6
Median	100.0
Min, Max	0, 100
Nights in hospital due to index lesion to 36 months	
N	1406
Mean ± SD	3.0 ± 7.8
Median	1.0
Min, Max	0, 165
Time to first clinically-driven TLR (days) through 60 months	
N	366
Mean ± SD	669.2 ± 462.5

Table 4. IN.PACT global/clinical cohort (5 years) (continued)

Parameter	IN.PACT global clinical cohort/IN.PACT DCB (N=1406 subjects) (N=3006 devices)
Median	571.5
Min, Max	1, 1751

^a Numbers are cumulative incidence % (number of failures) based on the Kaplan-Meier method.

Percutaneous transluminal angioplasty (PTA) for ISR – The primary clinical benefits include restoring lumen patency and blood flow and inhibition of restenosis. The measurable patient-relevant outcomes and quantitative clinical data from the Medtronic sponsored IN.PACT global study (DCB ISR cohort) at 1 year are shown in *Table 5*.

Table 5. IN.PACT global/DCB ISR cohort (1 year)

Parameter	IN.PACT Global/DCB ISR cohort (N=166 subjects) (N=386 devices)
Primary patency at 12 months	69.3% (79/114)
Primary sustained clinical improvement at 12 months	78.7% (118/150)
Clinically-driven TLR at 12 months	8.1% (13/160)
Acute success	
Device success (per device)	99.7% (385/386)
Procedural success (per lesion)	99.6% (231/232)
Clinical success (per subject)	98.2% (163/166)
Other major secondary endpoints	
Change in quality of life from baseline by EQ-5D index to 12 months	
N	140
Mean ± SD	0.1518 ± 0.3055
Median	0.1060
Min, Max	-0.741, 0.899
Walking impairment by WIQ (%) to 12 months	
N	145
Mean ± SD	75.3 ± 32.0
Median	100.0
Min, Max	0, 100
Nights in hospital due to index lesion to 12 months	
N	166
Mean ± SD	1.9 ± 3.9
Median	1.0
Min, Max	0, 48
Time to first clinically-driven TLR (days) through 12 months	
N	13
Mean ± SD	191.3 ± 122.0
Median	204.0
Min, Max	10, 350

PTA for obstructive lesions of the native AVF – The primary clinical benefits include restoring lumen patency and blood flow and inhibition of restenosis. The measurable patient-relevant outcomes and quantitative clinical data from the Medtronic sponsored IN.PACT AV Access study at 1 year are shown in *Table 6*.

Table 6. IN.PACT AV DCB (1 year)

Parameter	IN.PACT AV DCB (N=170 subjects) (N=212 devices)
Target lesion primary patency through 360 days — freedom from:	63.8% (90/141)
• CD-TLR	35.0% (49/140)
Access circuit thrombosis	2.9% (4/138)
Access circuit primary patency through 360 days — freedom from:	53.8% (78/145)
• Reintervention in access circuit	45.1% (65/144)
• Access circuit thrombosis	2.9% (4/138)
Acute success	
Device success (per device)	100.0% (212/212)
Procedural success (per subject)	73.5% (125/170)
Clinical success (per subject)	100.0% (159/159)

3.5 Intended users

IN.PACT Admiral DCB percutaneous transluminal angioplasty procedures should only be performed by physicians that have experience with interventional techniques in the vascular system. This device should only be used within a clinical setting in a sterilized environment.

3.6 Device performance characteristics

The IN.PACT Admiral DCB is a drug-coated, over-the-wire peripheral balloon catheter. The FreePac drug coating on the balloon consists of the drug paclitaxel and the excipient urea. Urea facilitates the release and transfer of paclitaxel into the arterial vessel wall. The balloon catheter physically dilates the vessel lumen by percutaneous transluminal angioplasty, and the drug coating is intended to reduce the proliferative response that is associated with restenosis. Paclitaxel stabilizes microtubules to reduce cell proliferation.

4 Warnings

- A signal for increased risk of late mortality has been identified following the use of paclitaxel-coated balloons and paclitaxel-eluting stents for femoropopliteal arterial disease beginning approximately 2-3 years post-treatment compared with the use of non-drug coated devices. There is uncertainty regarding the magnitude and mechanism for the increased late mortality risk, including the impact of repeat paclitaxel-coated device exposure. Physicians should discuss this late mortality signal and the benefits and risks of available treatment options with their patients.
- Inspect the IN.PACT Admiral DCB prior to the procedure to verify that the product is intact and functional. Do not use if the outer or the inner packaging is damaged or opened.
- Never apply positive pressure to the balloon during preparation.
- To reduce the potential for vessel damage, inflate the balloon to a diameter approximating that of the vessel just distal to the stenosis.

- Do not manipulate the inserted IN.PACT Admiral DCB without sufficient fluoroscopy.
- Do not withdraw the IN.PACT Admiral DCB from the lesion before completely deflating the balloon under vacuum.
- Do not expose the device to organic solvents, for example, alcohol.
- Do not manipulate the IN.PACT Admiral DCB while the balloon is inflated. The position of the balloon may only be changed with the guidewire in place.
- If resistance occurs during manipulation, ascertain the cause via fluoroscopy, road mapping, or digital subtraction angiography (DSA) before moving the IN.PACT Admiral DCB backward or forward.
- Do not under any circumstance move the guidewire during inflation of the IN.PACT Admiral DCB.
- Do not exceed the rated burst pressure (RBP). The RBP is based on the results of in vitro testing. At least 99.9% of IN.PACT Admiral DCB balloons (with 95% confidence) will not burst at or below their RBP. Use a pressure monitoring device to prevent over pressurization. Use of pressures higher than those specified on the product label may result in a ruptured balloon, causing possible intimal damage and dissection.
- Physicians using the IN.PACT Admiral DCB should have thorough PTA training and keep abreast of recent publications regarding PTA techniques.
- Do not use air or any gaseous medium to inflate the balloon. Use only the recommended inflation medium (equal parts contrast medium and saline solution).
- Do not use after the labeled Use By date.
- Do not use with Lipiodol™ or Ethiodol™ contrast media (or other such contrast media that incorporate components of these agents).

5 Precautions

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
 - Administer appropriate drug therapy (anticoagulant, vasodilator, etc.) to the patient according to standard protocols for PTA before inserting the IN.PACT Admiral DCB.
 - To minimize the introduction of air, aspirate and flush the system and keep a tight catheter connection throughout the procedure.
 - Take precautions to prevent or reduce clotting when any catheter is used. Flush or rinse all products entering the vascular system with sterile isotonic saline or a similar solution via the guidewire access port prior to use. Consider the use of systemic heparinization.
 - Use with caution for procedures involving calcified lesions due to the abrasive nature of these lesions.
 - Identify allergic reactions to contrast media, antiplatelet therapy, balloon catheters, and FreePac coating before treatment.
 - Catheter applications vary. Select the technique on the basis of the patient's condition and the experience of the interventionalist.
 - Never advance the IN.PACT Admiral DCB without the guidewire extending from the tip.
 - Store at controlled room temperature in a dry place. Keep away from sunlight.
 - Use a pressure-monitoring device to prevent overpressurization.
- Caution:** Larger sizes of the IN.PACT Admiral DCB may exhibit slower deflation times, particularly on long catheter shafts. For all balloon lengths except for 200 and 250 mm, allow approximately 60 seconds for full balloon deflation. For the 200 and 250 mm balloon lengths, allow approximately 120 seconds for full balloon deflation.
- After use, this product may be a biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable hospital, administrative, and government regulations.
 - The safety and effectiveness of using the IN.PACT Admiral DCB to treat peripheral vascular disease extending to or concomitant below the knee has not been established.

6 Potential complications/adverse effects

Possible adverse events which may be associated with the use of the IN.PACT Admiral DCB may include, but are not limited to:

- Abrupt vessel closure/thrombosis (acute total occlusion/reocclusion that may require surgical intervention)
- Access site pain, hematoma, hemorrhage, and/or local infection (bleeding may require transfusions)
- Allergic reaction to contrast medium, antiplatelet therapy, or catheter system components
- Aneurysm, pseudoaneurysm, or arteriovenous (AV) fistula
- Arrhythmias
- Balloon rupture
- Death
- Detachment of a component of the balloon and/or catheter system
- Dissection, perforation, or rupture of the artery
- Drug reactions
- Endocarditis
- Failure of the balloon to perform as intended (inflation/deflation/retrieval)
- Failure to deliver the balloon as intended (may release drug into unintended arterial segment)
- Hypotension/hypertension
- Ischemia/infarction of tissue/organ (severe ischemic events in treated limb may require amputation)
- Local or distal thromboembolic episodes
- Pain and tenderness at puncture sites
- Pyrogenic reaction
- Renal insufficiency or failure
- Restenosis of the dilated artery
- Sepsis/infection
- Short-term hemodynamic deterioration
- Systemic embolization
- Vessel spasms or recoil/prolonged arterial spasms

Potential adverse events not captured above, that may be unique to the paclitaxel drug coating include, but are not limited to:

- Allergic/immunologic reaction
- Alopecia
- Anemia
- Gastrointestinal symptoms
- Hematologic dyscrasias (including leucopenia, neutropenia, thrombocytopenia)
- Hepatic enzyme changes
- Histologic changes in vessel wall, including inflammation, cellular damage, or necrosis
- Myalgia/arthralgia
- Myelosuppression
- Peripheral neuropathy

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

7 Pharmacological interaction

Metabolic degradation of paclitaxel occurs in the liver via cytochrome P450 isoenzymes CYP2C8 and CYP3A4, resulting in the generation of 6-alpha-hydroxypaclitaxel and both 3'-p-hydroxypaclitaxel and 6-alpha, 3'-p-dihydroxypaclitaxel, respectively. Exercise caution when administering paclitaxel concomitantly with known substrates or inhibitors of CYP2C8 and CYP3A4 since no formal studies have been conducted with known competitors or inhibitors of these isoenzymes. Agents that could compete with or inhibit the activity of the CYP2C8 and CYP3A4 may increase paclitaxel plasma levels.

Formal drug interaction studies have not been conducted with the IN.PACT Admiral DCB. When deciding to use an IN.PACT Admiral DCB in a patient who is taking a drug with known interactions to paclitaxel, or when deciding to initiate therapy with such a drug in a patient who has recently been treated with an IN.PACT Admiral DCB, give consideration to the potential for both systemic and local drug interactions.

Please see Figure 3 (Chemical Structure of Paclitaxel) and Figure 4 (Chemical Structure of Urea) at the beginning of this Instructions for Use.

8 Summary of the paclitaxel meta-analysis and data from the IN.PACT DCB clinical studies

8.1 Summary of the paclitaxel meta-analysis

A meta-analysis of randomized controlled trials published in December 2018 by Katsanos et. al. identified an increased risk of late mortality at 2 years and beyond for paclitaxel-coated balloons and paclitaxel-eluting stents used to treat femoropopliteal arterial disease. In response to these data, the FDA performed a patient-level meta-analysis of long-term follow-up data from the pivotal premarket randomized trials of paclitaxel-coated devices used to treat femoropopliteal disease using available clinical data through May 2019. The meta-analysis also showed a late mortality signal in study subjects treated with paclitaxel-coated devices compared to patients treated with uncoated devices. Specifically, in the 3 randomized trials with a total of 1090 patients and available 5-year data, the crude mortality rate was 19.8% (range 15.9% to 23.4%) in patients treated with paclitaxel-coated devices compared to 12.7% (range 11.2% to 14.0%) in subjects treated with uncoated devices. The relative risk for increased mortality at 5 years was 1.57 (95% confidence interval 1.16 to 2.13), which corresponds to a 57% relative increase in mortality in patients treated with paclitaxel-coated devices. As presented at the June 2019 FDA Advisory Committee Meeting, an independent meta-analysis of similar patient-level data provided by VIVA Physicians, a vascular medicine organization, reported similar findings with a hazard ratio of 1.38 (95% confidence interval 1.06 to 1.80). Additional analyses have been conducted and are underway that are specifically designed to assess the relationship of mortality to paclitaxel-coated devices.

The presence and magnitude of the late mortality risk should be interpreted with caution because of multiple limitations in the available data, including wide confidence intervals due to a small sample size, pooling of studies of different paclitaxel-coated devices that were not intended to be combined, substantial amounts of missing study data, no clear evidence of a paclitaxel dose effect on mortality, and no identified pathophysiological mechanism for the late deaths.

Paclitaxel-coated balloons and stents improve blood flow to the legs and decrease the likelihood of repeat procedures to reopen blocked blood vessels compared to uncoated devices. The benefits of paclitaxel-coated devices (for example, reduced reinterventions) should be considered in individual patients along with potential risks (for example, late mortality).

8.2 Data from the IN.PACT DCB clinical studies

In the IN.PACT SFA IDE Trial, based on the analysis completed for the June 2019 FDA Advisory Committee Meeting using As-Treated cohort and vital status update, the Kaplan Meier cumulative mortality estimates at 2, 3, and 5 years are 7.3% [3.8%, 10.8%], 10.5% [6.4%, 14.6%], 15.7% [10.8%, 20.6%], respectively, for the IN.PACT Admiral DCB treatment device and 0.9% [0.2%, 2.7%], 2.8% [0.5%, 5.9%], 11.2% [5.3%, 17.1%], respectively, for the PTA control device.

A separate independent individual patient-level meta-analysis was performed and published by Schneider, et. al. in January 2019. The meta-analysis used data from patients treated with IN.PACT Admiral DCB. A total of 1,981 patients with up to 5-year follow-up were included in the meta-analysis. The purpose was to determine if a correlation exists between paclitaxel exposure and mortality. Data from 2 single-arm and 2 randomized independently adjudicated prospective studies of a paclitaxel DCB (n=1,837) and uncoated percutaneous transluminal angioplasty (PTA) (n=143) were included. Analyses of baseline, procedure, and follow-up data of individual patients were performed to explore correlations of paclitaxel dose with long-term mortality. Survival time by paclitaxel dose tercile was analyzed with adjustment of inverse probability weighting to correct baseline imbalances and study as random effect. A standard cohort compared DCB- and PTA-treated patients with similar characteristics by applying criteria from pivotal studies (n=712 DCB, n=143 PTA). A survival analysis stratified nominal paclitaxel dose by low, mid, and upper terciles; mean doses were 5,019.0, 10,007.5, and 19,978.2 µg, respectively. Rates of freedom from all-cause mortality between the 3 groups through 5 years were 85.8%, 84.2%, and 88.2%, respectively (p=0.731). There was no significant difference in all-cause mortality between DCB and PTA through 5 years comparing all patients (unadjusted p=0.092) or patients with similar characteristics (adjusted p=0.188). This independent patient-level meta-analysis demonstrates that the IN.PACT Admiral DCB is safe. Within DCB patients, there was no correlation between paclitaxel dose and mortality.

9 Instructions for use

Prior to angioplasty, carefully examine all equipment to be used during the procedure, including the IN.PACT Admiral DCB, to verify proper functioning. Verify that the IN.PACT Admiral DCB size is suitable for the specific procedure for which it is intended.

Handle the device with extreme caution in order to avoid any damage to the folded balloon. Avoid exposing the balloon's drug coating to excessive handling or contact with liquids prior to preparation and delivery as the coating may be susceptible to damage or premature drug release.

1. Inflation Device Preparation

- Prepare the inflation device according to the manufacturer's instructions.
- Prior to use, carefully examine the unit to verify that neither the IN.PACT Admiral DCB nor the sterile packaging have been damaged in shipment.

2. IN.PACT Admiral DCB Selection

- The nominal balloon size is equal to the inner diameter of the artery distal to the lesion. The balloon should extend approximately 1 cm beyond the lesion both proximally and distally. The use of several IN.PACT Admiral DCB devices in the same patient is possible in the case of long or multiple lesions. Additional balloons must overlap by 1 cm; do not overlap by more than 1 cm. If the stenosis cannot be crossed with the desired IN.PACT Admiral DCB, use a smaller catheter, without drug coating, to predilate the lesion and facilitate the passage of the appropriately sized IN.PACT Admiral DCB. In order to reduce potential procedure-related complications, use only the minimum number of devices needed to cover the lesions.

The safety and effectiveness of using multiple IN.PACT Admiral DCB devices with a total drug dosage up to 34,854 µg of paclitaxel has been evaluated in preclinical and clinical studies. Refer to Table 2 for nominal drug content by balloon size.

- In cases of total occlusion or subocclusive lesions, predilate with a non-drug-coated balloon.

3. IN.PACT Admiral DCB Preparation

- The IN.PACT Admiral DCB is packaged in a protective tray; carefully remove the catheter from the packaging.
- Purge the IN.PACT Admiral DCB of air bubbles prior to use. Keep the protective sheath in place during the purging procedure. Begin by pointing the IN.PACT Admiral DCB downward while holding it vertically. Connect a luer-lock syringe partially filled with the

saline-contrast mixture to the inflation port of the IN.PACT Admiral DCB hub. Apply negative pressure until air is completely evacuated and release the plunger. Repeat this operation until migration of air bubbles towards the syringe stops.

c. Remove the protective sheath from the balloon and discard. Do not use the protective sheath as an introduction aid or rewinding tool.

d. Flush the wire lumen properly through the luer-lock connector hub.

4. Inflation Device Connection to the IN.PACT Admiral DCB

- To remove any air lodged in the distal luer fitting of the inflation device, purge approximately 15 mL (cc) of contrast medium.
- With the stopcock in the closed position, disconnect the syringe used in preparation, slightly applying positive pressure. A meniscus of contrast medium will appear in the balloon port when the syringe is removed. Verify that a meniscus of contrast medium is evident in both the IN.PACT Admiral DCB balloon port hub and the inflation device's connection. Securely couple the inflation device to the balloon port of the IN.PACT Admiral DCB.

5. Use of the IN.PACT Admiral DCB

- Insert a guidewire through the hemostatic valve of the introducer sheath, following the manufacturer's instructions or standard practice. Advance the guidewire carefully into the introducer sheath.
 - Attach a torque device to the guidewire, if desired. Under fluoroscopy, advance the guidewire to the desired vessel, then across the stenosis. Remove the torque device once the guidewire is positioned.
 - Load the distal tip of the IN.PACT Admiral DCB onto the guidewire.
 - Open the hemostatic valve to allow for easy passage of the balloon and to prevent damage to the balloon coating. While the balloon is fully deflated, advance the IN.PACT Admiral DCB through the hemostatic valve of the introducer sheath. To avoid kinking, advance the IN.PACT Admiral DCB slowly, in small increments, until the proximal end of the guidewire emerges from the catheter. Close the hemostatic valve of the introducer sheath only as much as is needed to prevent blood return while permitting easy movements of the IN.PACT Admiral DCB. If the IN.PACT Admiral DCB encounters resistance, do not advance it through the adapter.
 - Under fluoroscopy, use the balloon's radiopaque markers to position the balloon within the lesion to be dilated, then inflate the balloon to the appropriate pressure (refer to balloon compliance table).
 - For optimal mechanical dilatation of the vessel, an inflation time of 180 seconds is strongly recommended for the IN.PACT Admiral DCB. Adequate drug transfer occurs in the first 60 seconds of inflation. In order to optimize lesion dilatation, longer inflation times are possible at the discretion of the operator.
- Note:** The IN.PACT Admiral DCB is intended for a single inflation only.
- Open the hemostatic valve and withdraw the completely deflated IN.PACT Admiral DCB through the hemostatic valve. Tighten the knurled knob on the hemostatic valve.
 - If necessary, the guidewire can remain in the vessel for the use of different balloon types or sizes.
 - WARNING: Do not dilate same vessel segment with more than one IN.PACT Admiral DCB. If the vessel requires postdilatation after treatment with the IN.PACT Admiral DCB, use a standard dilatation balloon without drug coating. Do not implant a drug-eluting stent into a vessel segment that has been treated by an IN.PACT Admiral DCB. Bare-metal stents are safe for use.**
 - Withdraw the guidewire from the introducer sheath through the hemostatic valve. When complete, withdraw the introducer sheath.
- Pre- and postprocedure antiplatelet regimen
 - Administer dual antiplatelet therapy (aspirin and clopidogrel or ticlopidine) before the procedure and for a minimum of 4 weeks after the intervention. Increase to 3 months if the patient is stented. Prolonged antiplatelet therapy can be given at the discretion of the physician.

10 Device disposal

Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

11 How supplied

The IN.PACT Admiral DCB is supplied sterile and for single use only. It is sterilized by ethylene oxide gas. If left unopened, this product will remain stable through the labeled Use By date. Do not use if the packaging is damaged.

12 Storage

Store the device in the original container. Store between 15 °C and 30 °C (59 °F and 86 °F). Use product by the Use-by Date noted on the package. Do not store near radiation or ultraviolet light sources.

13 Disclaimer of warranty

Note: This disclaimer of warranty does not apply in Australia.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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© 2024 Medtronic
Medtronic and Medtronic logo are trademarks of Medtronic. ^{TM*} Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic B.V.
Medtronic Limited