


STRUCTURAL HEART THERAPY

AMPLATZER™

SEPTAL OCCLUDERS



Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**.  
Always check the regulatory status of the device in your region.



THE PROVEN STANDARD

# FOR TRANSCATHETER ATRIAL SEPTAL DEFECT CLOSURE<sup>1, 2</sup>

Amplatzer™ Septal Occluders are the standard of care for minimally invasive atrial septal defect (ASD) closure. **These double-disc occluders are comprised of Nitinol mesh with polyester fabric.** They are designed to securely appose the septal wall on each side of the defect and create a platform for tissue in-growth after implantation.

## WHEN SIMPLICITY MATTERS

The primary treatment option for closure of secundum atrial septal defects is with transcatheter devices.<sup>2</sup> This is when trust and simplicity matter. The Amplatzer Septal Occluder is the proven standard of care in transcatheter ASD closure.<sup>1,2</sup>

### SAFETY IN NUMBERS

**98.5% Closure Rate:** With no significant residual shunt (less than 2 mm) at 1-year follow-up<sup>1</sup>

**Low Major and Minor Complication Rates:**

Postprocedure and long-term based on 5-year follow-up<sup>1</sup>


### SPECIFICALLY DESIGNED FOR ASD CLOSURE

**Wide Waist:** Centers device and fills the ASD<sup>3</sup>

**Shape-memory Nitinol Mesh:** Designed to securely appose both sides of the septal wall<sup>3</sup>

**Polyester Material:** Promotes occlusion and tissue in-growth<sup>3</sup>

**Precise Placement:** Device can be easily recaptured and redeployed<sup>3</sup>



The waist of the Amplatzer™ Septal Occluder fills the defect for optimal occlusion.

## WHEN FLEXIBILITY MATTERS

Not all septal defects are the same. They come in a variety of shapes, sizes and often there can be multiple communications between the left and right atria. These types of atrial septal defects require a special device designed with the flexibility to meet multiple needs.

The addition of the Amplatzer Multifenestrated Septal Occluder - “Cribriform” to the Amplatzer family of occluders enables transcatheter closure for the majority of atrial septal defects.


### SPECIFICALLY DESIGNED FOR MULTI-FENESTRATED ASD CLOSURE

**Narrow Waist:** Allows for placement through a central defect<sup>3</sup>

**Matched Disc Diameters:** Maximizes coverage of multiple fenestrations<sup>3</sup>

**Polyester Material:** Promotes occlusion and tissue in-growth<sup>3</sup>

**Precise Placement:** Device can be easily recaptured and redeployed<sup>3</sup>



The Amplatzer™ Cribriform Multi-Fenestrated Septal Occluder enables occlusion of the defect by covering the fenestrations with a single device.

Stocking both **Amplatzer™ Septal and Cribriform Occluders** provides you the ability to choose the appropriate device for a wide size-range of defects.

# ORDERING INFORMATION

## AMPLATZER™ SEPTAL OCCLUDER

Model / Reorder Number	Waist Diameter (mm)	Waist Width (mm)	Right Atrial Disc Diameter (mm)	Left Atrial Disc Diameter (mm)
9-ASD-004	4	3	12	16
9-ASD-005	5	3	13	17
9-ASD-006	6	3	14	18
9-ASD-007	7	3	15	19
9-ASD-008	8	3	16	20
9-ASD-009	9	3	17	21
9-ASD-010	10	3	18	22
9-ASD-011	11	4	21	25
9-ASD-012	12	4	22	26
9-ASD-013	13	4	23	27
9-ASD-014	14	4	24	28
9-ASD-015	15	4	25	29
9-ASD-016	16	4	26	30
9-ASD-017	17	4	27	31
9-ASD-018	18	4	28	32
9-ASD-019	19	4	29	33
9-ASD-020	20	4	30	34
9-ASD-022	22	4	32	36
9-ASD-024	24	4	34	38
9-ASD-026	26	4	36	40
9-ASD-028	28	4	38	42
9-ASD-030	30	4	40	44
9-ASD-032	32	4	42	46
9-ASD-034	34	4	44	50
9-ASD-036	36	4	46	52
9-ASD-038	38	4	48	54
9-ASD-040	40	4	50	56

## AMPLATZER™ CRIBRIFORM MULTI-FENESTRATED SEPTAL OCCLUDER

Model / Reorder Number	Right & Left Atrium Disc Diameter (mm)	Waist Width (mm)
9-ASD-MF-018	18	3
9-ASD-MF-025	25	3
9-ASD-MF-030	30	3
9-ASD-MF-035	35	3
9-ASD-MF-040	40	3

### REFERENCE:

1. Amplatzer Septal Occluder Instructions for Use.
2. Kashour TS, Latroche B, Elhoury ME, et al. Successful Percutaneous Closure of a Secundum Atrial Septal Defect through Femoral Approach in a Patient with Interrupted Inferior Vena Cava. Congenital Heart Disease. 2010;5(6):620-623.
3. Test(s) performed by and data on file at Abbott.

**CAUTION:** This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at [eifu.abbottvascular.com](http://eifu.abbottvascular.com) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Always check the regulatory status of the device in your region.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott.

### Abbott

3200 Lakeside Dr, Santa Clara, CA 95054, USA  
[www.cardiovascular.abbott](http://www.cardiovascular.abbott)

™ Indicates a trademark of the Abbott group of companies.

© 2020 Abbott. All Rights Reserved. MAT-2004327 v1.0 | Item approved for Global OUS use only.






**ST. JUDE MEDICAL™**

M8007064 Rev. F

Declaration of Conformity

## SJM Declaration of Conformity

St. Jude Medical (SJM) hereby declares that the following SJM facilities and products conform to the applicable provisions of Annex II of the Medical Device Directive (MDD) 93/42/EEC. All supporting documentation is retained under the premises of SJM. We declare no application has been lodged with any other notified body for the same products. This declaration is issued under the sole responsibility of the manufacturer. This declaration supersedes any declaration issued previously for the same product(s).

**Manufacturer Address:**

AGA Medical Corporation  
5050 Nathan Lane North  
Plymouth, Minnesota 55442, USA

**European Representative:**

St. Jude Medical Coordination Center BVBA  
The Corporate Village  
Da Vincilaan 11 Box F1  
1935 Zaventem, Belgium

**Product Type:**

Cardiac Occluder

**Product Name(s):**

AMPLATZER Septal Occluder  
AMPLATZER Multi-Fenestrated Septal Occluder -  
"Cribiform"  
AMPLATZER PFO Occluder

**Model Number(s):**

Product	Model Numbers			Original CE Mark Date
AMPLATZER Septal Occluder	9-ASD-004	9-ASD-013	9-ASD-024	24 February 1998
	9-ASD-005	9-ASD-014	9-ASD-026	
	9-ASD-006	9-ASD-015	9-ASD-028	
	9-ASD-007	9-ASD-016	9-ASD-030	
	9-ASD-008	9-ASD-017	9-ASD-032	
	9-ASD-009	9-ASD-018	9-ASD-034	
	9-ASD-010	9-ASD-019	9-ASD-036	
	9-ASD-011	9-ASD-020	9-ASD-038	
	9-ASD-012	9-ASD-022	9-ASD-040	
AMPLATZER Multi-Fenestrated Septal Occluder – "Cribiform"	9-ASD-MF-018	9-ASD-MF-035		2 September 2002
	9-ASD-MF-025	9-ASD-MF-040		
	9-ASD-MF-030			
AMPLATZER PFO Occluder	9-PFO-018	9-PFO-025		24 February 1998
	9-PFO-030	9-PFO-035		

**Signature:**

  
Lisa Becker  
Senior Director, Regulatory Affairs

**Issue Date**

  
2 Feb 18



## SJM Declaration of Conformity

<b>Classification:</b>	Class III per Annex II, Rule 8
<b>GMDN Code(s):</b>	45418
<b>EC Design Certificate No and Expiration Date:</b>	Certificate No: CE 594291 Expiration Date: 23 February 2023
<b>Annex II, Clause 3 Certificate No and Expiration Date:</b>	Certificate No: CE 590631 Expiration Date: 23 February 2023
<b>Applicable Quality System Standards:</b>	ISO 13485
<b>Notified Body:</b>	BSI Kitemark Court Davy Avenue Knowlhill Milton Keynes MK5 8PP UK
<b>Notified Body Number:</b>	0086
<b>AMPLATZER Septal Occluder Manufacturing Facilities:</b>	AGA Medical Corporation 5050 Nathan Lane North Plymouth, Minnesota 55442 USA  St. Jude Medical, Costa Rica Ltda. Edificio #44, Calle 0, Avenida 2 Zona Franca Coyol, El Coyol, Alajuela, Costa Rica
<b>AMPLATZER Multi-Fenestrated Septal Occluder - "Cribiform" Manufacturing Facility:</b>	AGA Medical Corporation 5050 Nathan Lane North Plymouth, Minnesota 55442 USA
<b>AMPLATZER PFO Occluder Manufacturing Facility:</b>	AGA Medical Corporation 5050 Nathan Lane North Plymouth, Minnesota 55442 USA

**Signature:**

  
Lisa Becker  
Senior Director, Regulatory Affairs



Issue Date

# EC Design-Examination Certificate

Directive 93/42/EEC on Medical Devices, Annex II Section 4

**No.****CE 594291****Issued To:**

**AGA Medical Corporation  
5050 Nathan Lane North  
Plymouth  
Minnesota  
55442  
USA**

**In respect of:**

**AMPLATZER Septal Occluder, AMPLATZER Multifenestrated Septal Occluder - "Cribiform" and  
AMPLATZER PFO Occluder**

BSI has performed a design examination on the above devices in accordance with the Council Directive 93/42/EEC, Annex II Section 4. The design conforms to the requirements of this directive. For marketing of these products an additional Annex II excluding Section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 0086):



Stewart Brain, Head of Compliance & Risk -  
Medical Devices

First Issued: **2013-02-24**

Date: **2018-02-14**

Expiry Date: **2023-02-23**

...making excellence a habit.™

# EC Design-Examination Certificate

## Supplementary Information to CE 594291

Issued To:

**AGA Medical Corporation**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

### Amplatzer Septal Occluder:

Device Size (mm)	Model Number	Wire Diameter (mm)/Number of Wires	Distal Disc Diameter (mm)	Proximal Disc Diameter (mm)	Waist Diameter	Waist Length
4	9-ASD-004	.004/72	16	12	4	3
5	9-ASD-005	.004/72	17	13	5	3
6	9-ASD-006	.004/72	18	14	6	3
7	9-ASD-007	.004/72	19	15	7	3
8	9-ASD-008	.004/72	20	16	8	3
9	9-ASD-009	.004/72	21	17	9	3
10	9-ASD-010	.004/72	22	18	10	3
11	9-ASD-011	.005/72	25	21	11	4
12	9-ASD-012	.005/72	26	22	12	4
13	9-ASD-013	.005/72	27	23	13	4
14	9-ASD-014	.005/72	28	24	14	4
15	9-ASD-015	.005/72	29	25	15	4
16	9-ASD-016	.005/72	30	26	16	4

First Issued: **2013-02-24**

Date: **2018-02-14**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 2 of 5

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
 This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP. Tel: + 44 345 080 9000  
 BSI Assurance UK Limited, registered in England under number 7805321 at 389 Chiswick High Road, London W4 4AL, UK.  
 A member of BSI Group of Companies.

# EC Design-Examination Certificate

## Supplementary Information to CE 594291

Issued To:

**AGA Medical Corporation**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

Device Size (mm)	Model Number	Wire Diameter (mm)/Number of Wires	Distal Disc Diameter (mm)	Proximal Disc Diameter (mm)	Waist Diameter	Waist Length
17	9-ASD-017	.005/72	31	27	17	4
18	9-ASD-018	.006/72	32	28	18	4
19	9-ASD-019	.006/72	33	29	19	4
20	9-ASD-020	.006/72	34	30	20	4
22	9-ASD-022	.006/72	36	32	22	4
24	9-ASD-024	.006/72	38	34	24	4
26	9-ASD-026	.007/72	40	36	26	4
28	9-ASD-028	.007/72	42	38	28	4
30	9-ASD-030	.007/72	44	40	30	4
32	9-ASD-032	.007/72	46	42	32	4
34	9-ASD-034	.008/72	50	44	34	4
36	9-ASD-036	.008/72	52	46	36	4
38	9-ASD-038	.008/72	54	48	38	4
40	9-ASD-040	.008/72	56	50	40	4

First Issued: **2013-02-24**

Date: **2018-02-14**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 3 of 5

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
 This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP. Tel: + 44 345 080 9000  
 BSI Assurance UK Limited, registered in England under number 7805321 at 389 Chiswick High Road, London W4 4AL, UK.  
 A member of BSI Group of Companies.

# EC Design-Examination Certificate

## Supplementary Information to CE 594291

Issued To:

**AGA Medical Corporation**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

### Amplatzer Multi-Fenestrated Septal Occluder – “Cribriform”:

Device Size (mm)	Model Number	Wire Diameter (mm)/Number of Wires	Disc Diameter (mm)	Diameter of Raised Disc (mm)
18	9-ASD-MF-018	.005/72	18	6
25	9-ASD-MF-025	.005/72	25	13
30	9-ASD-MF-030	.005/72	30	30
35	9-ASD-MF-035	.006/72	35	23
40	9-ASD-MF-040	.006/72	40	40

### AMPLATZER PFO Occluder:

Model Number	Device Size (mm)	Right Atrial Disc Diameter (mm)	Left Atrial Disc Diameter (mm)	Waist Length
9-PFO-018	18	18	18	3
9-PFO-025	25	25	18	3
9-PFO-030	30	30	30	3
9-PFO-035	35	35	25	3

First Issued: **2013-02-24**

Date: **2018-02-14**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 4 of 5

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
 This certificate was issued electronically and is bound by the conditions of the contract.



# EC Design-Examination Certificate

## Supplementary Information to CE 594291

Issued To:

**AGA Medical Corporation**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

## Certificate History

Date	Reference Number	Action
24 February 2013	10139223	First Issue – Transfer from another Notified Body.
28 May 2013	10141409	Change the surface finish of the nitinol wire from black oxide to chemically etched.
03 February 2014	10144456	Update to Amplatzer Septal Occluder IFU. There are no changes to the intended use.
06 February 2015	10152724	Addition of St. Jude Medical Costa Rica Ltda. as an alternate manufacturing site for the Amplatzer Septal Occluder. Addition of Synergy Health AST, SRL as an alternate sterilization site for the Amplatzer Septal Occluder.
25 August 2015	10156926	DuPont Tyvek Medical Transition Project update.
01 February 2016	10160623	Addition of Sterigenics Willowbrook, IL as a sterilizer.
Current	8872315	Certificate renewal.

First Issued: **2013-02-24**

Date: **2018-02-14**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 5 of 5

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP. Tel: + 44 345 080 9000  
BSI Assurance UK Limited, registered in England under number 7805321 at 389 Chiswick High Road, London W4 4AL, UK.  
A member of BSI Group of Companies.



AMPLATZER™ TREVISIO™  
INTRAVASCULAR DELIVERY SYSTEM

MORE FLEXIBILITY,  
MORE CONTROL



DISCOVER THE UNIQUE  
AMPLATZER™ TREVISIO™

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.





# RELIABLE PRECISION WHEN IT MATTERS MOST

The Amplatzer™ Trevisio™ Intravascular Delivery System is an ultra-flexible delivery system enabling interventional cardiologists to perform their work with complete confidence. It leverages the one-piece cable design utilized by the Amplatzer™ TorqVue™ Delivery System, also known as the Classic Amplatzer™ Delivery System\*. Trevisio is designed for no compromises on torque strength, sheath diameter and pushability.



## ULTRA-FLEXIBLE TIP

The diagram illustrates the Amplatzer Trevisio delivery system. A long, thin, curved catheter is shown with a yellow line tracing its path. Two circular callouts provide detailed views of specific components. The top callout shows the ultra-flexible tip, which has a series of rounded, overlapping segments. The bottom callout shows the flexible transition section, which has a series of rounded, overlapping segments. The main catheter body is shown with a series of rounded, overlapping segments, and the tip is shown with a series of rounded, overlapping segments.

- Improves assessment of device position before cable release.
- Reduces bias on the device.

## FLEXIBLE TRANSITION SECTION

Maintains sheath position during deployment of the device.

## CLASSIC AMPLATZER™ DELIVERY SYSTEM\*



## AMPLATZER™ TREVISIO™ INTRAVASCULAR DELIVERY SYSTEM


ULTRA-  
FLEXIBLE TIP

FLEXIBLE  
TRANSITION  
SECTION

STIFF  
PROXIMAL  
SECTION



Flexible tip reduces the bias on the device and improves the assessment of device position prior to cable release.



### STIFF PROXIMAL SECTION

Maintains pushability of the  
delivery system.

## DELIVERY SYSTEM DIMENSIONS

DELIVERY SYSTEM (SHEATH SIZE)	INNER DIAMETER OF SHEATH	OUTER DIAMETER OF SHEATH	MODEL NUMBER/DELIVERY SYSTEM SIZE (MM)
6 Fr	2.11 mm (0.08 in)	2.79 mm (0.11 in)	9-ATV06F45/60
7 Fr	2.44 mm (0.10 in)	3.18 mm (0.13 in)	9-ATV07F45/60
7 Fr	2.44 mm (0.10 in)	3.18 mm (0.13 in)	9-ATV07F45/80
8 Fr	2.69 mm (0.11 in)	3.45 mm (0.14 in)	9-ATV08F45/60
8 Fr	2.69 mm (0.11 in)	3.45 mm (0.14 in)	9-ATV08F45/80
9 Fr	3.00 mm (0.12 in)	3.81 mm (0.15 in)	9-ATV09F45/80
10 Fr	3.30 mm (0.13 in)	4.14 mm (0.16 in)	9-ATV10F45/80
12 Fr	3.99 mm (0.16 in)	4.80 mm (0.19 in)	9-ATV12F45/80
13 Fr	4.32 mm (0.17 in)	5.13 mm (0.20 in)	9-ATV13F45/80

\*Amplatzer™ TorqVue™ Delivery System



## DOWNLOAD NOW

Download the free Amplatzer™ Portfolio App from your favorite app store.



**CAUTION:** This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at [eifu.abbottvascular.com](http://eifu.abbottvascular.com) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs.

### Abbott

Park Lake, Culliganlaan 2b, 1831 Diegem, Belgium  
[www.cardiovascular.abbott](http://www.cardiovascular.abbott)

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third party trademark, which is property of its respective owner.

© 2020 Abbott. All Rights Reserved. MAT-2001046 v1.0 | Item approved for OUS use only.



# AMPLATZER™ TREVISIO™ INTRAVASCULAR DELIVERY SYSTEM


## DEVICE DESCRIPTION

The Amplatzer™ Trevisio™ Intravascular Delivery System is an ultra-flexible delivery system enabling interventional cardiologists to perform their work with complete confidence. It leverages the one-piece cable design utilized by the Amplatzer™ TorqVue™ Delivery System, also known as the Classic Amplatzer™ Delivery System. Trevisio is designed for no compromises on torque strength, sheath diameter and pushability.

The Amplatzer™ Trevisio™ Intravascular Delivery System comes with a 45° sheath curve in 60 or 80 cm lengths from 6Fr to 13Fr.

## INDICATIONS AND USAGE


The Amplatzer™ Trevisio™ Intravascular Delivery System is intended to facilitate the attachment, loading, delivery and deployment of the Amplatzer™ Occluder devices. The Amplatzer™ Trevisio™ Intravascular Delivery System is indicated for use by interventional cardiologists to place an Amplatzer™ Occluder device.



## DEVICE SPECIFICATIONS

Model Number/ Delivery System Size mm	Delivery System (Sheath Size) Fr	Inner Diameter of Sheath mm (inch)	Outer Diameter of Sheath mm (inch)	Length cm
9-ATV06F45/60	6	2.11 (0.08)	2.79 (0.11)	60
9-ATV07F45/60	7	2.44 (0.10)	3.18 (0.13)	60
9-ATV07F45/80	7	2.44 (0.10)	3.18 (0.13)	80
9-ATV08F45/60	8	2.69 (0.11)	3.45 (0.14)	60
9-ATV08F45/80	8	2.69 (0.11)	3.45 (0.14)	80
9-ATV09F45/80	9	3.00 (0.12)	3.81 (0.15)	80
9-ATV10F45/80	10	3.30 (0.13)	4.14 (0.16)	80
9-ATV12F45/80	12	3.99 (0.16)	4.80 (0.19)	80
9-ATV13F45/80	13	4.32 (0.17)	5.13 (0.20)	80

## PACKAGING INFORMATION



- A. Loader – Introduces an Amplatzer™ device into the sheath
- B. Hemostasis valve with extension tube and stopcock – Allows flushing of the delivery system and controls blood backflow
- C. Sheath – Provides a pathway through which an Amplatzer™ device is delivered
- D. Dilator – Eases penetration of tissue and minimizes vessel trauma
- E. Delivery cable – Attaches to the device to control its movement through the sheath
- F. Plastic vise – Attaches to the delivery cable and serves as a handle for disconnecting (unscrewing) the delivery cable from a device

Note that the device and the delivery system are being sold separately.

## ASSOCIATED AMPLATZER™ PRODUCTS

Amplatzer™ Trevisio™ Intravascular Delivery System Sizes							
	6 Fr	7 Fr	8 Fr	9 Fr	10 Fr	12 Fr	13 Fr
Amplatzer™ Septal (ASD) Occluder	9-ASD-004	9-ASD-011	9-ASD-018	9-ASD-020	9-ASD-026	9-ASD-032	9-ASD-032
	9-ASD-005	9-ASD-012	9-ASD-019	9-ASD-022	9-ASD-028	9-ASD-034	9-ASD-034
	9-ASD-006	9-ASD-013	–	9-ASD-024	9-ASD-030	9-ASD-036	9-ASD-036
	9-ASD-007	9-ASD-014	–	–	–	9-ASD-038	9-ASD-038
	9-ASD-008	9-ASD-015	–	–	–	9-ASD-040	9-ASD-040
	9-ASD-009	9-ASD-016	–	–	–	–	–
	9-ASD-010	9-ASD-017	–	–	–	–	–
Amplatzer™ Multi-fenestrated ASD (Cribriform) Occluder	–	–	9-ASD-MF-018			–	–
	–	–	9-ASD-MF-025	9-ASD-MF-035	9-ASD-MF-040	–	–
	–	–	9-ASD-MF-030			–	–
Amplatzer™ PFO Occluder	–	–	9-PFO-018		–	–	–
	–	–	9-PFO-025	9-PFO-035	–	–	–
	–	–	9-PFO-030		–	–	–
Amplatzer™ Muscular VSD Occluder	9-VSD-MUSC-004				–	–	–
	9-VSD-MUSC-006	9-VSD-MUSC-012	9-VSD-MUSC-014	9-VSD-MUSC-018	–	–	–
	9-VSD-MUSC-008		9-VSD-MUSC-016		–	–	–
	9-VSD-MUSC-010				–	–	–
Amplatzer™ Muscular PI VSD Occluder	–	–	–	9-VSDMUSCPI-016	9-VSDMUSCPI-020	–	–
	–	–	–	9-VSDMUSCPI-018	9-VSDMUSCPI-022	–	–
	–	–	–		9-VSDMUSCPI-024	–	–

## LATEX-FREE INFORMATION

These Amplatzer™ products do not contain latex.

### Abbott Vascular International BVBA

Park Lane, Culliganlaan 2B, B-1831 Diegem, Belgium, Tel: +32 2 714 14 11

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at [eifu.abbottvascular.com](http://eifu.abbottvascular.com) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott. Information contained herein for DISTRIBUTION in Europe, Middle East and Africa ONLY. Please check the regulatory status of the device before distribution in areas where CE marking is not the regulation in force.

For more information, visit our web site at [www.abbott.com](http://www.abbott.com)

© 2020 Abbott. All rights reserved. 9-EH-11038-01 03-2020



## Abbott Medical Declaration of Conformity

Abbott Medical (Abbott) hereby declares that the following Abbott facilities and products conform to the applicable provisions of Annex II of the Medical Device Directive (MDD) 93/42/EEC. All supporting documentation is retained under the premises of Abbott. We declare no application has been lodged with any other notified body for the same products. This declaration is issued under the sole responsibility of the manufacturer. This declaration supersedes any declaration issued previously for the same product(s).

**Manufacturer Address:**

Abbott Medical  
5050 Nathan Lane North  
Plymouth, Minnesota 55442, USA

**European Representative:**

Abbott Medical  
The Corporate Village  
Da Vincilaan 11 Box F1  
1935 Zaventem, Belgium

**Product Type:**

Cardiac Occluder Delivery Kit

**Product Name(s):**

Amplatzer TorqVue Delivery System  
Amplatzer TorqVue Exchange System  
Amplatzer TorqVue 2 Delivery Sheath  
Amplatzer TorqVue LP Delivery System  
Amplatzer TorqVue LP Catheter  
Amplatzer TorqVue Delivery System with Pusher Catheter  
Amplatzer TorqVue 45°x45° Delivery Sheath  
Amplatzer Amulet Delivery Sheath  
Amplatzer Trevisio Intravascular Delivery System

**Signature:**

Michelle Grossman  
Director, Regulatory Affairs

April 21, 2020  
Issue Date




## Abbott Medical Declaration of Conformity

Model Number(s):

Product	Model Numbers	Original CE Mark Date
Amplatzer TorqVue Delivery System	9-ITV06F45/60      9-ITV13F45/80 9-ITV07F45/60      9-ITV05F180/60 9-ITV07F45/80      9-ITV06F180/60 9-ITV08F45/60      9-ITV06F180/80 9-ITV08F45/80      9-ITV07F180/80 9-ITV09F45/80      9-ITV08F180/80 9-ITV10F45/80      9-ITV09F180/80 9-ITV12F45/80	10 October 2005
Amplatzer TorqVue Exchange System	9-EITV09F45/80 9-EITV12F45/80 9-EITV06F180/80 9-EITV08F180/80 9-EITV10F180/80	10 October 2005
Amplatzer TorqVue 2 Delivery Sheath	9-TV2-05F120 9-TV2-06F120 9-TV2-07F120	19 February 2010
Amplatzer TorqVue LP Delivery System	9-TVLP4F90/060 9-TVLP4F90/080 9-TVLP5F90/060 9-TVLP5F90/080	07 February 2008
Amplatzer TorqVue LP Catheter	9-TVLPC4F90/080	28 April 2011
Amplatzer TorqVue Delivery System with Pusher Catheter	9-ITVP07F-180/80 9-ITVP08F-180/80 9-ITVP09F-180/80	21 June 2011
Amplatzer TorqVue 45°x45° Delivery Sheath	9-TV45X45-09F-100 9-TV45X45-10F-100 9-TV45X45-12F-100 9-TV45X45-13F-100 9-TV45X45-14F-100	03 December 2008 (9-13 Fr) 24 February 2012 (14 Fr)
Amplatzer Amulet Delivery Sheath	DS-TV45X45-12F-080 DS-TV45X45-14F-080	08 February 2017
Amplatzer Trevisio Intravascular Delivery System	9-ATV06F45/60      9-ATV09F45/80 9-ATV07F45/60      9-ATV10F45/80 9-ATV07F45/80      9-ATV12F45/80 9-ATV08F45/60      9-ATV13F45/80 9-ATV08F45/80	20 April 2020

Signature:

  
 \_\_\_\_\_  
 Michelle Grossman  
 Director, Regulatory Affairs

*April 21 2020*

 \_\_\_\_\_  
 Issue Date

## Abbott Medical Declaration of Conformity

<b>Classification:</b>	Class III (Rule 7) Annex II, Section 4 GHTF Class D
<b>GMDN Code(s):</b>	45419
<b>EC Design Certificate No and Expiration Date:</b>	Certificate No: CE 694956 Expiration Date: 23 February 2023
<b>Annex II, Clause 3 Certificate No and Expiration Date:</b>	Certificate No: CE 694788 Expiration Date: 23 February 2023
<b>Applicable Quality System Standards:</b>	ISO 13485
<b>Notified Body:</b>	BSI Group The Netherlands B.V. Say Building John M. Keynesplein 9 1066 EP Amsterdam The Netherlands
<b>Notified Body Number:</b>	2727

**Signature:**  
\_\_\_\_\_  
Michelle Grossman  
Director, Regulatory AffairsApril 21, 2020  
Issue Date

# EC Design-Examination Certificate

Directive 93/42/EEC on Medical Devices, Annex II Section 4

**No.** **CE 694956**  
**Issued To:** **Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

In respect of:

**Amplatzer Delivery Systems**

BSI has performed a design examination on the above devices in accordance with the Council Directive 93/42/EEC, Annex II Section 4. The design conforms to the requirements of this directive. For marketing of these products an additional Annex II excluding Section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):



Gary E Slack, Senior Vice President Medical Devices

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 1 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
 BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
 A member of BSI Group of Companies.

# EC Design-Examination Certificate

## Supplementary Information to CE 694956

Issued To:

**Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

<b>Amplatzer™ TorqVue™ Delivery System</b>			
<b>Intended purpose per IFU:</b>			
The Amplatzer™ TorqVue™ Delivery System is intended to facilitate the attachment, loading, delivery and deployment of the Amplatzer™ Occluder devices.			
<b>Classification:</b> Class III			
Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-ITV06F45/60	6	45°	60
9-ITV07F45/60	7	45°	60
9-ITV07F45/80	7	45°	80
9-ITV08F45/60	8	45°	60
9-ITV08F45/80	8	45°	80
9-ITV09F45/80	9	45°	80
9-ITV10F45/80	10	45°	80
9-ITV12F45/80	12	45°	80
9-ITV13F45/80	13	45°	80
9-ITV05F180/60	5	180°	60
9-ITV06F180/60	6	180°	60
9-ITV06F180/80	6	180°	80
9-ITV07F180/80	7	180°	80
9-ITV08F180/80	8	180°	80
9-ITV09F180/80	9	180°	80

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 2 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
A member of BSI Group of Companies.

# EC Design-Examination Certificate

## Supplementary Information to CE 694956

Issued To:

**Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

<b>Amplatzer™ TorqVue™ Exchange System</b>			
<b>Intended purpose per IFU:</b> The Amplatzer™ TorqVue™ Exchange System is intended for removal of an Amplatzer™ Delivery Sheath and subsequent exchange for an Amplatzer™ Delivery Sheath of equal or larger diameter.			
<b>Classification:</b> Class III			
Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-EITV09F45/80	9	45°	80
9-EITV12F45/80	12	45°	80
9-EITV06F180/80	6	180°	80
9-EITV08F180/80	8	180°	80
9-EITV10F180/80	10	180°	80

<b>Amplatzer™ TorqVue™ 2 Delivery Sheath</b>			
<b>Intended purpose per IFU:</b> The Amplatzer™ TorqVue™ 2 Delivery Sheath is intended to provide a pathway through which devices are introduced within the chambers and coronary vasculature of the heart or in the peripheral vasculature.			
<b>Classification:</b> Class III			
Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-TV2-05F120	5	none	120
9-TV2-06F120	6	none	120
9-TV2-07F120	7	none	120

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 3 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
 This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
 BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
 A member of BSI Group of Companies.



# EC Design-Examination Certificate

## Supplementary Information to CE 694956

Issued To:

**Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

<b>Amplatzer™ TorqVue™ Delivery System with Pusher Catheter</b>			
<b>Intended purpose per IFU:</b> The Amplatzer™ TorqVue™ Delivery System with Pusher Catheter is intended to facilitate the attachment, loading, delivery and deployment of the Amplatzer™ Membranous VSD Occluder device.			
<b>Classification:</b> Class III			
Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-ITVP07F180/80	7	180°	80
9-ITVP08F180/80	8	180°	80
9-ITVP09F180/80	9	180°	80

<b>Amplatzer™ TorqVue™ LP Delivery System</b>				
<b>Intended purpose per IFU:</b> The Amplatzer™ TorqVue™ LP Delivery System is intended to facilitate the attachment, loading, delivery, and deployment of the Amplatzer™ devices.				
<b>Classification:</b> Class III				
Catalogue Number	Model, Type			
	Device Size (Fr)	Curve Dimension	Length (cm)	Delivery Wire Length (cm)
9-TVLP4F90/060	4	90°	60	160
9-TVLP4F90/080	4	90°	80	195
9-TVLP5F90/060	5	90°	60	160
9-TVLP5F90/080	5	90°	80	195

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 4 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
A member of BSI Group of Companies.

# EC Design-Examination Certificate

## Supplementary Information to CE 694956

Issued To:

**Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

### Amplatzer™ Trevisio™ Intravascular Delivery System

#### Intended purpose per IFU:

The Amplatzer™ Trevisio™ Intravascular Delivery System is intended to facilitate the attachment, loading, delivery and deployment of the Amplatzer™ Occluder devices.

**Classification:** Class III

Catalogue Number	Model, Type		
	Sheath Size (Fr)	Tip Angle	Usable Length (cm)
9-ATV06F45/60	6	45°	60
9-ATV07F45/60	7	45°	60
9-ATV07F45/80	7	45°	80
9-ATV08F45/60	8	45°	60
9-ATV08F45/80	8	45°	80
9-ATV09F45/80	9	45°	80
9-ATV10F45/80	10	45°	80
9-ATV12F45/80	12	45°	80
9-ATV13F45/80	13	45°	80

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 5 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780

BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.

A member of BSI Group of Companies.



# EC Design-Examination Certificate

## Supplementary Information to CE 694956

Issued To:

**Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

<b>Amplatzer™ TorqVue™ LP Catheter</b>				
<b>Intended purpose per IFU:</b> The TorqVue™ LP Catheter is intended to facilitate the loading, delivery, and deployment of Amplatzer™ devices.				
<b>Classification:</b> Class III				
Catalogue Number	Model, Type			
	Device Size (Fr)	Usable Length (cm)	Tip Outer Diameter mm (in)	Tip Inner Diameter mm (in)
9-TVLPC4F90/080	4	80	1.40 (.055)	1.17 (0.046)

<b>Amplatzer™ TorqVue™ 45x45 Delivery Sheath</b>		
<b>Intended purpose per IFU:</b> The Amplatzer™ TorqVue™ Delivery Sheath is intended to provide a pathway through which devices are introduced within the chambers of the heart.		
<b>Classification:</b> Class III		
Catalogue Number	Model, Type	
	Sheath Size (Fr)	Length (cm)
9-TV45X45-09F-100	9	100
9-TV45X45-10F-100	10	100
9-TV45X45-12F-100	12	100
9-TV45X45-13F-100	13	100
9-TV45X45-14F-100	14	100

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 6 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
 This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
 BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
 A member of BSI Group of Companies.

# EC Design-Examination Certificate

## Supplementary Information to CE 694956

Issued To:

**Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

### Amplatzer™ Amulet™ Delivery Sheath

#### Intended purpose per IFU:

The Amplatzer™ Amulet™ Delivery Sheath is intended to provide a pathway through which devices are introduced within the chambers of the heart.

**Classification:** Class III

Catalogue Number	Model, Type	
	Sheath Size (Fr)	Length (cm)
DS-TV45X45-12F-080	12	80
DS-TV45X45-14F-080	14	80

First Issued: **2018-09-03**

Date: **2020-04-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 7 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
A member of BSI Group of Companies.

# EC Design-Examination Certificate

## Supplementary Information to CE 694956

Issued To:

**Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

## Certificate History

Date	Reference Number	Action
03 September 2018	8957249	First Issue. Mirror certificate to CE 594294.
20 February 2019	8243107	Traceable to NB 0086.
20 March 2019	9738457	Addition of Sterigenics US LLC, Salt Lake City, Utah for ETO Sterilization.
16 December 2019	3053900	Addition of Midwest Sterilization Corporation, Jackson, Missouri USA for ETO Sterilization in chambers 1, 2, 3, 6, and 13.
Current	9784335	Extension to scope to include the Amplatzer™ Trevisio™ Intravascular Delivery System. Revision to scope statement to remove "TorqVue" brand name. Reformat product tables. Correction of delivery wire length for Amplatzer™ TorqVue™ LP Delivery System. Correction of tip outer diameter for Amplatzer™ TorqVue™ LP Catheter.

First Issued: **2018-09-03**

Date: **2020-04-20**


Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 8 of 8

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
A member of BSI Group of Companies.



## AMPLATZER™ Sizing Balloon II

Accessory

### Product Highlights

- Triple lumen balloon catheter designed to measure cardiovascular structures
- Compliant balloon material supports low pressure inflation and allows precise measurement of a variety of defect shapes when using stop-flow technique
- Three radiopaque marker bands strategically placed to assist with alignment of imaging and support proper sizing calibration radiographically and echocardiographically
- Soft distal tip and flexible shaft offer smooth access for atraumatic balloon positioning
- Ultra thick membrane requires no dilatation of the entry site; balloon can be introduced over the guidewire

### Ordering Information

Contents: 1 sizing balloon

Model/Reorder Number	Maximum Defect Size	Maximum Inflation Volume	Balloon Length	Shaft Size	Usable Length	Guidewire
9-SB-018	20 mm	12 cc	3.5 cm	6 F	70 cm	0.035 inch
9-SB-024	27 mm	25 cc	4.5 cm	7 F	70 cm	0.035 inch
9-SB-034	40 mm	90 cc	5.5 cm	8 F	70 cm	0.035 inch

## AMPLATZER™ Sizing Plate

Accessory

### Ordering Information

Contents: 1 sizing plate

Model/Reorder Number	Description
9-ASD-SZP	Measuring plate with circular openings ranging from 4 to 38 mm



## SJM Declaration of Conformity

St. Jude Medical (SJM) hereby declares that the following SJM facilities and products conform to the applicable provisions of Annex II of the Medical Device Directive (MDD) 93/42/EEC. All supporting documentation is retained under the premises of SJM. We declare no application has been lodged with any other notified body for the same products. This declaration is issued under the sole responsibility of the manufacturer. This declaration supersedes any declaration issued previously for the same product(s).

<b>Manufacturer Address:</b>	AGA Medical Corporation 5050 Nathan Lane North Plymouth, Minnesota 55442, USA
<b>European Representative:</b>	St. Jude Medical Coordination Center BVBA The Corporate Village Da Vincilaan 11 Box F1 1935 Zaventem, Belgium
<b>Product Type:</b>	Cardiac Catheter, Balloon, Sizing
<b>Product Name(s):</b>	AMPLATZER Sizing Balloon II
<b>Model Number(s):</b>	9-SB-018, 9-SB-024, 9-SB-034
<b>Classification:</b>	Class III (Rule 7) Annex II, Section 4, GHTF Class D
<b>GMDN Code(s):</b>	42417
<b>Original CE Mark Date:</b>	10 June 2005 (24 mm) 02 Sept 2005 (18, 34 mm)
<b>EC Design Certificate No and Expiration Date:</b>	Certificate No: CE 595439 Expiration Date: 23 Feb 2023
<b>Annex II, Clause 3 Certificate No and Expiration Date:</b>	Certificate No: CE 590631 Expiration Date: 23 Feb 2023
<b>Applicable Quality System Standards:</b>	ISO 13485
<b>Notified Body:</b>	BSI Kitemark Court Knowlhill Milton Keynes MK5 8PP UK
<b>Notified Body Number:</b>	0086

Signature:

  
Lisa Becker  
Senior Director, Regulatory Affairs

  
Issue Date



By Royal Charter

# EC Design-Examination Certificate

Directive 93/42/EEC on Medical Devices, Annex II Section 4

**No.****CE 595439****Issued To:****AGA Medical Corporation  
5050 Nathan Lane North  
Plymouth  
Minnesota  
55442  
USA**

In respect of:

**AMPLATZER Sizing Balloon II**

BSI has performed a design examination on the above devices in accordance with the Council Directive 93/42/EEC, Annex II Section 4. The design conforms to the requirements of this directive. For marketing of these products an additional Annex II excluding Section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 0086):

Stewart Brain, Head of Compliance & Risk -  
Medical Devices

First Issued: **2013-02-24**Date: **2018-01-10**Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 1 of 3





By Royal Charter

# EC Design-Examination Certificate

## Supplementary Information to CE 595439

Issued To:

**AGA Medical Corporation**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

Finished Product Number	Balloon Diameter (mm)	Volume at Diameter (cc)	Useable Balloon Length (mm)	Catheter Useable Length (cm)	Shaft French Size (F)	Maximum Recommended Inflation Volume (cc)	Minimum Volume at Failure (cc)
9-SB-018	18	10	35	70	6	12	25
9-SB-024	24	25	45	70	7	25	60
9-SB-034	34	55	55	70	8	90	120

First Issued: **2013-02-24**Date: **2018-01-10**Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 2 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP. Tel: + 44 345 080 9000





By Royal Charter

# EC Design-Examination Certificate

## Supplementary Information to CE 595439

Issued To:

**AGA Medical Corporation**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

## Certificate History

Date	Reference Number	Action
24 February 2013	10139846	First Issue – Transfer from another Notified Body.
25 August 2015	10156926	DuPont Tyvek Medical Transition Project update.
01 February 2016	10160623	Addition of Sterigenics Willowbrook, IL as a sterilizer.
Current	8864699	Certificate renewal.

First Issued: **2013-02-24**Date: **2018-01-10**Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 3 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.

This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Kitemark Court, Davy Avenue, Knowlhill, Milton Keynes MK5 8PP. Tel: + 44 345 080 9000

# AMPLATZER™ GUIDEWIRES

## Ordering Information

Contents: 1 guidewire

Model/ Reorder No	Diam (in)	Body	Floppy Tip Lgth (cms)	Tip Description	Usable Lgth (cm)
9-GW-001	0.035	Super Stiff	6	7.5 mm, Modified J-tip	260
9-GW-002	0.035	Super Stiff	5	1.5 mm, Modified J-tip	260
9-GW-003	0.035	Super Stiff	20	6 mm, J-tip	300
9-GW-004*	0.035	Soft Tip, Fixed Core	NA	6 mm, J-tip	300

\*9-GW-004, also referred to as “Noodlewire,” is a soft tip, flexible guidewire recommended for establishing an arterial-venous loop, facilitating closure of ventricular septal defects.



# Amplatzer™ Guidewire

**Guidewire**

es: Guia  
pt: Fio-Guia  
tr: Kılavuz Teli

UDI



(01)00811806010724(17)261231(10)9941803



**2026-12-31**


REF

9-GW-002



LOT

9941803



	0.035/0.889
	260 cm
	1.5 mm MOD J
	PTFE
	FIXED



STERILE EO



Does not contain  
natural rubber  
latex components



2022-01-24

CE  
2797



**Abbott Medical**  
5050 Nathan Lane North  
Plymouth, MN  
55442 USA  
+1 855 478 5833  
+1 651 756 5833

EC REP

**Abbott Medical**  
The Corporate Village  
Da Vincilaan 11 Box F1  
1935 Zaventem  
Belgium  
+32 2 774 68 11

**Manufactured By**

Lake Region Medical Limited,  
Butlersland, New Ross, Co.  
Wexford, Ireland.

Made in Ireland



AGABF-031V

600035212 B

**SJM Declaration of Conformity**

St. Jude Medical (SJM) hereby declares that the following SJM facilities and products conform to the applicable provisions of Annex II of the Medical Device Directive (MDD) 93/42/ECC. All supporting documentation is retained under the premises of SJM. We declare no application has been lodged with any other notified body for the same products. This declaration is issued under the sole responsibility of the manufacturer. This declaration supersedes any declaration issued previously for the same product(s).

**Manufacturer Address:** AGA Medical Corporation  
5050 Nathan Lane North  
Plymouth, Minnesota 55442, USA

**European Representative:** St. Jude Medical Coordination Center BVBA  
The Corporate Village  
Da Vincilaan 11 Box F1  
1935 Zaventem, Belgium

**Product Type:** Catheter Guide Wire

**Product Name(s):** AMPLATZER Guidewires

**Model Number(s):** 9-GW-001, 9-GW-002, 9-GW-003, 9-GW-004

**Classification:** Class III (Rule 7) Annex II, Section 4  
GHTF Class D

**GMDN Code(s):** 35094

**Original CE Mark Date:** 22 March 2001 (1-3), 22 Sept 2003 (4, Noodlewire)

**EC Certificate No and expiration date:** Certificate No: CE 594293  
Expiration Date: 23 Feb 2023


**Annex II, Clause 3:** Certificate No: CE 590631  
Expiration Date: 23 Feb 2023

**Applicable Quality System Standards:** ISO 13485

**Notified Body:** BSI  
Kitemark Court  
Davy Avenue  
Knowlhill  
Milton Keynes  
MK5 8PP  
UK

**Notified Body Number:** 0086

**Signature:**

  
Lisa Becker  
Senior Director, Regulatory Affairs

  
Issue Date




# EC Design-Examination Certificate

Directive 93/42/EEC on Medical Devices, Annex II Section 4

**No.****CE 694955****Issued To:****Abbott Medical  
5050 Nathan Lane North  
Plymouth  
Minnesota  
55442  
USA****In respect of:****AMPLATZER Guidewires**

BSI has performed a design examination on the above devices in accordance with the Council Directive 93/42/EEC, Annex II Section 4. The design conforms to the requirements of this directive. For marketing of these products an additional Annex II excluding Section 4 certificate is required.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 2797):



Albert Roossien, Regulatory Lead

**First Issued: 2018-09-03****Date: 2019-02-20****Expiry Date: 2023-02-23**

...making excellence a habit.™

Page 1 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
A member of BSI Group of Companies.



# EC Design-Examination Certificate

## Supplementary Information to CE 694955

Issued To:

**Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

Part Number	Description
9-GW-001	JFC-SS Modified J, Fixed, PTFE Coated, Super Stiff
9-GW-002	J1.5FC-SS, Modified J, Fixed, PTFE Coated, Super Stiff
9-GW-003	J9FC-FS-LLLT Fixed Core, Long (20 cm) PTFE Coated, Finger-Straightenable
9-GW-004	Noodlewire Guidewire

First Issued: **2018-09-03**

Date: **2019-02-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 2 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
A member of BSI Group of Companies.

# EC Design-Examination Certificate

## Supplementary Information to CE 694955

Issued To:

**Abbott Medical**  
**5050 Nathan Lane North**  
**Plymouth**  
**Minnesota**  
**55442**  
**USA**

## Certificate History

Date	Reference Number	Action
03 September 2018	8957249	First Issue. Mirror certificate to CE 594293.
Current	8243107	Traceable to NB 0086.

First Issued: **2018-09-03**

Date: **2019-02-20**

Expiry Date: **2023-02-23**

...making excellence a habit.™

Page 3 of 3

Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive as demonstrated through the required surveillance activities of the Notified Body.  
This certificate was issued electronically and is bound by the conditions of the contract.

Information and Contact: BSI, Say Building, John M. Keynesplein 9, 1066 EP Amsterdam, The Netherlands Tel: + 31 20 346 0780  
BSI Group The Netherlands B.V. registered in The Netherlands under 33264284.  
A member of BSI Group of Companies.