Total T4

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Total T4

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of total circulating thyroxine (T4) in serum. Measurements of T4 are used in the diagnosis and treatment of thyroid disease.

Catalog Numbers:

L2KT42 (200 tests), **L2KT46** (600 tests)

Test Code: **TT4** Color: **Light Blue**

Summary and Explanation

Thyroxine (T4), the principal thyroid hormone, normally circulates at a level of approximately 4.5 to 12.5 µg/dL (58 to 161 nmol/L), largely bound to transport proteins, especially TBG. Given normal levels of thyroid hormone-binding proteins, hyperthyroidism is characterized by increased levels of circulating T4, hypothyroidism by decreased levels.

Exceptions to this parallelism between thyroid status and total T4 concentration are found, however, in patients with abnormal levels of thyroid-binding proteins. Levels of TBG are known to be altered under various physiological, pharmacological and genetic conditions. Thus, elevated T4 levels may be obtained when TBG levels are high, as in pregnancy, acute intermittent porphyria, hyperproteinemia, hereditary TBG elevation and in patients undergoing estrogen therapy or taking oral contraceptives. Total T4 levels may be depressed when TBG levels are low, as in nephrotic, hepatic, gastrointestinal and neoplastic disorders; in acromegaly, hypoproteinemia and hereditary TBG deficiency; and in patients undergoing androgen, testosterone or anabolic steroid therapy. Diphenylhydantoin and large doses of salicylates and liothyronine may also cause low T4 values (not reflective of thyroid status) due to their competition for binding sites on TBG.

Because the total T4 level is often abnormal in euthyroid subjects and because it may be normal in subjects with

abnormal thyroid function, an estimate of the level of circulating TBG is desirable, as given for example by a T3 Uptake assay. In disorders of thyroid function, the total T4 and the T3 Uptake values will deviate from normal in the same direction, whereas in euthyroid patients with the TBG alterations they will deviate in opposite directions. The product of the T4 and the T3 Uptake values, divided by 100, is known as the Free Thyroxine Index (FT4I), a widely used indicator of thyroid status.

Principle of the Procedure

IMMULITE 2000 Total T4 is a solid-phase, competitive chemiluminescent enzyme immunoassay.

Incubation Cycles: 1 × 30 minutes

Specimen Collection

EDTA plasma should not be used as a sample type.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Total T4 has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 15 µL serum

Storage: 7 days at 2–8°C, or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.¹³⁻¹⁵

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Total T4 Bead Pack (L2T412)

With barcode. 200 beads, coated with monoclonal murine anti-T4. Stable at 2–8°C until expiration date.

L2KT42: 1 pack **L2KT46:** 3 packs

Total T4 Reagent Wedge (L2T4A2)

With barcode. 21 mL alkaline phosphatase (bovine calf intestine) conjugated to T4, with preservative, dispensed equally into chambers A and B. Stable at 2–8°C until expiration date.

L2KT42: 1 wedge **L2KT46:** 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Total T4 Adjustors (LT4L, LT4H)

Two vials (Low and High), 2 mL each, of T4 in processed human serum, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KT42: 1 set **L2KT46:** 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

Multi-Diluent 1 (L2M1Z)

For the on-board dilution of high samples. One vial of concentrated (ready-to-use) processed, normal human serum, containing undetectable to low levels of total T4, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2M1Z: 25 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M1Z: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

Also Required

Distilled or deionized water; test tubes; controls.

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of T4.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Expected Values

Adult: Based on its relationship to IMMULITE Total T4 (see Method Comparison), the assay can be expected to have essentially the same reference ranges for adults:

4.5–12.5 µg/dL (58–161 nmol/L)

Reference ranges were generated using serum samples collected in glass tubes without gel barriers or clot-promoting additives.

Pediatric: Reference intervals for the pediatric population (children and adolescents) were established for the IMMULITE Total T4 assay in accordance with CLSI guideline EP28-A3C.¹⁶

For analysis of data, the population was divided into three age subgroups:

- Infants: subjects aged 1–23 months
- Children: subjects aged 2–12 years
- Adolescents: subjects aged 13–20 years

A non-parametric approach was used to establish the reference intervals for children and adolescents where the 2.5 and 97.5 percentiles of the distribution of values were calculated. For the infant population, a robust measure of location and spread, as developed by Horn and Pesce, was used to estimate the 2.5 and the 97.5 percentile reference intervals, accommodating the smaller sample size.^{16–18}

The reference intervals detailed by age group and number of samples are presented in the Reference Intervals table.

IMMULITE 2000 Total T4 Pediatric Reference Intervals

Age Group	<i>n</i>	Conventional (µg/dL)	SI (nmol/L)
Infants (1–23 Months)	82	6.2 – 11.8	80 – 152
Children (2–12 Years)	197	5.4 – 11.1	70 – 143
Adolescents (13–20 Years)	148	4.9 – 10.2	63 – 131

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

Many physiological, pharmacological, pathological, and genetic factors affect the interpretation of total T4 results. (Refer to the Summary and Explanation section.)

More than 20% of patients with non-thyroidal disease who are critically ill have low serum T4 levels.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in

combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in µg/dL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

µg/dL × 12.87 → nmol/L

Calibration Range: 1.0 to 24 µg/dL
(13 to 309 nmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 0.3 µg/dL
(3.9 nmol/L)

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Linearity: Samples were assayed under various dilutions. (See "Linearity" table for representative data.)

Recovery: Samples spiked 1 to 19 with three T4 solutions (42, 90 and 150 µg/dL), were assayed. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for T4. (See "Specificity" table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 512 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 5000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 20 volunteers into plain, heparinized and Becton Dickinson SST® vacutainer tubes. Equal volumes of the matched samples were spiked with various concentrations of Total T4, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 Total T4 procedure.

(Heparin) = 0.89 (Serum) + 0.61 µg/dL
r = 0.987

(SST) = 1.01 (Plain Tubes) + 0.42 µg/dL
r = 0.971

Means:

7.72 µg/dL (Serum)

7.51 µg/dL (Heparin)

8.21 µg/dL (SST)

Method Comparison: The assay was compared to IMMULITE Total T4 on 276 samples. (Concentration range: approximately 3 to 24 µg/dL. See graph.)
By linear regression:

(IML 2000) = 0.93 (IML) + 0.56 µg/dL
r = 0.978

Means:

9.2 µg/dL (IMMULITE 2000)

9.3 µg/dL (IMMULITE)

References

- 1) Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982.
- 2) Britton KE, Quinn V, Brown BL, Ekins RP. A strategy for thyroid function tests. Br Med J 1975; iii:350–2.
- 3) Ekins RP, editor. Free thyroid hormones. Amsterdam: Excerpta Medica, 1979.
- 4) Gershengorn MC, Glinoe D, Robbins J. Transport and metabolism of thyroid hormones. In: Visscher M De, editor. The thyroid gland. New York: Raven Press, 1980: 81–121.
- 5) Prince HP, Ramsden DB. A new theoretical description of the binding of thyroid hormones by serum proteins. Clin Endocrinol 1977;7:307–24.
- 6) Refetoff S. Thyroid hormone transport. In: DeGroot LJ, editor. Endocrinology. Philadelphia: Grune and Stratton, 1979:1:347–56.
- 7) Refetoff S. Thyroid function tests. In: DeGroot LJ, editor. Endocrinology, Philadelphia: Grune and Stratton, 1979:1:387–428.
- 8) Robbins J, Johnson ml. Theoretical considerations in the transport of thyroid hormones in blood. In: Ekins RP, editor. Free thyroid hormones. Amsterdam: Excerpta Medica, 1979: 1–16.
- 9) Robbins J, Rall JE. The iodine-containing hormones. In: Gray CH, James VHT, editors. Hormones in blood. London: Academic Press, 1979: 575–688.
- 10) Slag MF, Morley JE, Elson MK, et al. Hypothyroxinemia in critically ill patients as a predictor of high mortality. JAMA 1981;245: 43–5.
- 11) Walfish PG. Thyroid physiology and

pathology. In: Collu R Jr, Ducharme H, Guyda, editors. Pediatric endocrinology. New York: Raven Press, 1981: 357–31. 12)Wenzel KW. Pharmacological interference with *in vitro* tests of thyroid function. Metabolism 1981;30:717–32. 13) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377–82, 387–8. 14) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3. 15) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030. 16) Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Guideline EP28-A3C. 17) Horn PS, Pesce AJ, Reference Intervals. A User's Guide, Washington, DC:AACC Press; 2005. 18) Reed AH, Henry RJ, Mason WB. Influence of statistical method used on the resulting estimate of normal range. 1971; 17:275-284.

Technical Assistance

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Tables and Graphs

Precision (µg/dL)

	Mean ³	Within-Run ¹			Total ²	
		SD ⁴	CV ⁵	SD	CV	
1	1.8	0.20	11.1%	0.21	11.7%	
2	2.6	0.23	8.9%	0.28	10.8%	
3	5.2	0.32	6.2%	0.44	8.5%	
4	7.0	0.37	5.3%	0.43	6.1%	
5	8.2	0.38	4.6%	0.46	5.6%	
6	13	0.60	4.6%	0.78	6.0%	
7	16	0.74	4.6%	0.89	5.6%	

Recovery (µg/dL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	2.7	—	—
	A	5.1	4.7	109%
	B	7.8	7.1	110%
	C	8.5	10.1	84%
2	—	5.5	—	—
	A	8.3	7.3	114%
	B	10.2	9.7	105%
	C	13	13	100%
3	—	7.8	—	—
	A	11	9.5	116%
	B	12	12	100%
	C	14	15	93%

Linearity (µg/dL)


	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	8 in 8 ⁵	8.9	—	—
	4 in 8	4.2	4.5	93%
	2 in 8	2.0	2.3	87%
	1 in 8	1.1	1.1	100%
2	8 in 8	9.6	—	—
	4 in 8	4.9	4.8	102%
	2 in 8	2.7	2.4	113%
	1 in 8	1.4	1.2	117%
3	8 in 8	10.6	—	—
	4 in 8	5.6	5.3	106%
	2 in 8	2.7	2.7	100%
	1 in 8	1.5	1.3	115%
4	8 in 8	11.4	—	—
	4 in 8	5.7	5.7	100%
	2 in 8	2.7	2.9	93%
	1 in 8	1.5	1.4	107%

Specificity (µg/dL)

Compound ¹	Amount Added ²	% Cross reactivity ³
L-Thyroxine (T4)	—	100%
D-Thyroxine	10	74%
Tetraiodothyroacetic acid	10	50%
Triiodo-L-thyronine (T3)	25	1.7%
Reverse T3	25	ND
Triiodo-D-thyronine	1	ND
Triiodothyroacetic acid	500 10	1.6% ND
Monoiodotyrosine	1000	ND
Diiodo-L-tyrosine	1000	ND
Methimazole	1000	ND
5,5'-Diphenylhydantoin	1000	ND
Phenylbutazone	1000	ND
6- <i>n</i> -Propyl-2-thiouracil	1000	ND

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2018 Siemens Healthcare Diagnostics. All rights reserved.
CONTROL AB