HbA1c

mindray

Hemoglobin A1c Kit (Enzymatic Assay Method)

Order Information

Cat. No.	Package size
105-002165-00	R(Hb): 1×30 mL + R1(HbA1c): 1×30 mL + R2(HbA1c): 1×12 mL + R0: 1×150 mL
105-002166-00	R(Hb): 1×40 mL + R1(HbA1c): 1×40 mL + R2(HbA1c): 1×15 mL + R0: 1×200 mL
105-002167-00	R(Hb): 2×40 mL + R1(HbA1c): 2×40 mL + R2(HbA1c): 2×15 mL + R0: 2×150 mL
105-009338-00	R1: 1×40 mL + R2: 1×15 mL + R0: 3×40 mL
105-005737-00	R1: 1×40 mL + R2: 1×15 mL + R0: 1×200 mL
105-005738-00	R1: 2×40 mL + R2: 2×15 mL + R0: 3×58 mL

Intended Purpose

In vitro test for the quantitative determination of Hemoglobin A1c (HbA1c) concentration in human whole blood on Mindray BS series chemistry analyzers. It is intended to be used for diagnosis and monitoring therapeutic effect of diabetes.


Summary¹⁻⁷

Hemoglobin (Hb) consists of four protein chains with four heme portions, and is the red-pigmented protein located in the erythrocytes. Its main function is the transport of oxygen and carbon dioxide in blood. Each Hb molecule is able to bind four oxygen molecules. Hb consists of a variety of subfractions and derivatives. Among this heterogeneous group of hemoglobins HbA1c is one of the glycosylated hemoglobins, a subfraction formed by the attachment of various sugars to the Hb molecule. HbA1c is formed in two steps by the non-enzymatic reaction of glucose with the N-terminal amino group of the β -chain of normal adult Hb (HbA). The first step is reversible and yields labile HbA1c. This slowly rearranges in the second reaction step to yield stable HbA1c. In the erythrocytes, the relative amount of HbA converted to stable HbA1c increases with the average concentration of glucose in the blood. The conversion to stable HbA1c is limited by the erythrocyte's life span of approximately 100 to 120 days. As a result, HbA1c reflects the average blood glucose level during the preceding 2 to 3 months. HbA1c is thus suitable to monitor long-term blood glucose control in individuals with diabetes mellitus. More recent glucose levels have a greater influence on the HbA1c level. The approximate relationship between HbA1c and mean blood glucose value during the preceding 2 to 3 months has been analyzed by several studies.

HbA1c

Assay Principle

Enzymatic Assay Method

In the first reaction, the concentration of hemoglobin is measured at absorbance of fixed wavelength, and simultaneously the fructosyl dipeptides are generated from the N-terminus amino groups of the beta-chain of HbA1c by the reaction of proteinase. In the second reaction, the reaction of Fructosyl peptide oxidase (FPOX) with fructosyl dipeptides, generated hydroperoxide allows 10-(carboxymethylaminocarbonyl)-3,7-bis (dimethyl-amino) phenothiazine sodium salt (DA67) to develop a color in the presence of peroxidase. The change in absorbance is measured for HbA1c determination. The combined assay results for hemoglobin and HbA1c are used by the system to calculate and express HbA1c%.

Reagents Components

	N,N-Dimethylformamide	<10.0 g/L
	Calcium chloride	<2.0 g/L
R(Hb)/R1(HbA1c)/R1:	Proteinase	> 500 KU/L
	DA67	0.1 mmol/L
	ProClin 300	0.05%
	Sodium azide	<0.01%
R2(HbA1c)/R2:	Peroxidase	>1500 U/L
	Fructosyl peptide oxidase	>1500 U/L
R0 (Pretreatment Solution)	NaNO ₂	1.0 g/L
	ProClin 300	0.05%

Storage and stability

Up to expiration date indicated on the label, when stored unopened at 2-8°C and protected from light.

On board in use, the reagents are stable for 28 days when refrigerated on the analyzer.

Contamination must be avoided.

Do not freeze the reagent.

HbA1c

These sections are applied for HbA1c assayed using whole blood samples with Centrifugation.


Specimen collection and preparation

■ Specimen types

K₂-EDTA, K₃-EDTA, NaF-EDTA, lithium heparin and sodium heparin whole blood are suitable for samples.

■ Preparation for Analysis

1. Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.
2. Prior to testing, whole blood samples should be centrifuged at 2000 rpm for 5 min.
3. Refer to the flowchart below:

■ On-board lysis

1. Refer to "BS-series System Operations Manual" -> "General Operating Procedure" -> "Whole Blood Test".
2. The proportion of centrifuged blood corpuscle sample and Pretreatment Solution is 10ul:200ul.

- Warning: For the anticoagulation tubes, the height of total sample in the tube should be no higher than 55mm and the height of centrifuged blood corpuscle should be no lower than 10mm. If not satisfied, adjust to the height of total sample and centrifuged blood corpuscle or change to the manual lysis.
- Warning: Micro-sample cups and Eppendorf microfuge tube are not allowed.

■ Manual lysis

- Get 25 μ L blood corpuscle from the centrifuged deposition into a sample cup or an Eppendorf microfuge tube. The blood corpuscle is easily adsorbed to the sample probe and pipette, which may cause exceed the target volume.
- Add 500 μ L of Pretreatment Solution into test tube.
- Close the test tube and lysis the blood by inverting 10 times or shaking for 10s on the vortex mixer.
- The hemolysate can be used as working samples after 4 minutes.
- The volume of blood corpuscle and Pretreatment Solution can be amplified in equal proportion if necessary.
- Specimens should be tested as soon as possible after pre-analytical treatment.

■ Sample Stability

- Whole blood Stability: 8 hours at 15-25°C
7 days at 2-8°C^B
30 days at (-25)-(-15)°C
- Hemolysate Stability: 8 hours at 15-25°C
24 hours at 2-8°C

Do not freeze the whole blood samples. Sample stability claims were established by manufacturer and/or based upon references, each laboratory should establish its own sample stability criteria.

Reagent Preparation

R(Hb)/R1(HbA1c)/R1, R2(HbA1c)/R2 and Pretreatment Solution are ready to use.

Please perform scheduled maintenance and standard operation including calibration and analysis to assure the performance of measurement system.

Materials required but not provided

- General laboratory materials: distilled/deionized water.
- Calibrator and Control: Please check the section of reagent instruction of Calibration and Quality Control.
- Mindray BS series chemistry analyzers and General laboratory equipment.

HbA1c

Assay procedure

■ Non-twin test

The Hb was tested by using reagent R(Hb), and the HbA1c was tested by using reagent R1 (HbA1c)/R2(HbA1c).

■ Twin test

The Hb and HbA1c were simultaneously tested by using reagent R1/R2.

1.Hb

Parameters Item	BS-800 chemistry analyzers
Assay type	Endpoint
Wavelength (Primary/Secondary)	505/800 nm
Reaction direction	Increase
R (Hb)/R1(HbA1c)	180 μ L
incubate at 37°C for 1 min, read the absorbance A1, then add:	
Sample or Calibrator	12 μ L
Mix thoroughly, incubate at 37°C for 5 min, read the absorbance A2, Then calculate $\Delta A_{\text{Hb-calibrator}}=(A2-A1)$ calibrator, $\Delta A_{\text{Hb-sample}}=(A2-A1)$ sample	

2.HbA1c

Parameters Item	BS-800 chemistry analyzers
Assay type	Endpoint
Wavelength (Primary/Secondary)	660/800 nm
Reaction direction	Increase
R1 (HbA1c)/R1	180 μ L
Sample or Calibrator	12 μ L
Mix, incubate at 37°C for 5 min, read the absorbance A1, then add:	
R2 (HbA1c)/R2	60 μ L
Mix thoroughly, incubate at 37°C for 5 min, read the absorbance A2, Then calculate $\Delta A_{\text{HbA1c-calibrator}}=(A2-A1)$ calibrator, $\Delta A_{\text{HbA1c-sample}}=(A2 - A1)$ sample	

3.HbA1c%

HbA1c% is the most common form in the reports. Report results of hemoglobin A1c in mmol/mol (IFCC) or percent (NGSP):

Protocol 1 according to IFCC: $\text{HbA1c (mmol/mol)} = \text{HbA1c/Hb} \times 1000$

Protocol 2 according to NGSP/DCCT: $\text{HbA1c(\%)} = 91.5 \times \text{HbA1c / Hb} + 2.15$

Parameters may vary in different chemistry analyzers, may adjust in proportion if necessary. For Mindray BS series chemistry analyzers, Reagent Parameters is available on request. Please refer to the appropriate operation manual for the analyzers.

Calibration

1. It is recommended to use the Mindray Calibrator (HbA1c Calibrator: 105-003680-00 or other suitable calibrators) for two-point calibration. Traceability of the HbA1c Calibrator can refer to the calibrator instructions for use of Mindray Company.
2. Calibration frequency
Calibration is stable for approximately 28 days on BS-800 chemistry analyzers. The calibration stability may vary on different instruments, each laboratory should set a calibration frequency in the instrument parameters appropriate to their usage pattern.
Recalibration may be necessary when the following occur:
 - As changed reagent lot.
 - As required following quality control procedures or out of control.
 - As executes specific maintenance or troubleshooting procedure of chemistry analyzers.
3. The calibrator values are lot-specific with the matched models listed in the value sheet.

Quality control

1. It is recommended to use the Mindray Control (HbA1c Control: 105-002140-00, 105-002138-00 or other suitable controls) to verify the performance of the measurement procedure; other suitable control material can be used in addition.
2. Two levels of control material are recommended to analyze each batch of samples. In addition, the control should be run with each new calibration, with each new reagent cartridge, and after specific maintenance or trouble shooting procedures as detailed in the appropriate system manual.
3. Each laboratory should establish its own internal quality control scheme and procedures for corrective action if control doesn't recover within the acceptable tolerances.

Calculation

$C_{\text{sample}} = (\Delta A_{\text{sample}} / \Delta A_{\text{calibrator}}) \times C_{\text{calibrator}}$

The BS series chemistry analyzer detects the change of absorbance(ΔA) and calculates the Hb/HbA1c concentration of each sample automatically after calibration.

The system calculates HbA1c and Hb results in $\mu\text{mol/L}$ (SI units). To convert the results to g/dL (common units), use the following equation:

Conversion factor: $\mu\text{mol/L} \times 0.1354 = \text{g/dL}$.

Expected values⁹

	Sample Type	S.I. Units
Whole blood	According to IFCC	20.20-42.06 mmol/mol
	According to NGSP/DCCT	4.0%-6.0%

The expected value is provided from reference, and Mindray has verified it by 60 samples of people from southern China.

Each laboratory should establish its own reference intervals based upon its particular locale and population characteristics since expected values may vary with geography, race, sex and age.

Performance Characteristics

■ Analytical Sensitivity

The Hemoglobin A1c Kit has an analytical sensitivity of 60 $\mu\text{mol/L}$ for Hb, 2 $\mu\text{mol/L}$ for HbA1c on BS-400. Analytical sensitivity is defined as the lowest concentration of analyte that can be distinguished from a sample that contains no analyte. It is calculated as the value lying 3 standard deviations above that of the mean from 20 replicates of an analyte-free sample.

■ Measuring range

The Mindray BS series systems provide the following linearity range:

Sample Type	S.I. Units
Whole blood	3%-16% (NGSP/DCCT)

A high concentration HbA1c sample (approximately 16%) is mixed with a low concentration sample (<3%) at different ratios, generating a series of dilutions. The HbA1c concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient $r \geq 0.9900$. The reportable range is 3%-16% (NGSP/DCCT), as the ratio can't be solved through dilution.

■ Precision

Precision was determined by following CLSI Approved Guideline EP05-A3¹⁰, each sample was assayed 2 times per run, 2 runs per day, a total 20 of days. The precision data (calculated according to NGSP/DCCT) of controls on BS-400 are summarized in the following tables*.

Specimen Type (N=80)	Mean (%)	Repeatability		Within-Lab	
		SD (%)	CV %	SD (%)	CV %
Control N	4.43	0.14	3.19	0.17	3.79
Control P	9.32	0.13	1.43	0.16	1.75

**Representative data, results in different instruments, laboratories may vary.*

■ Analytical Specificity

The samples with different concentration interfering substance were prepared by addition of interferent to human sample pools, and recovers are within ± 10.0 % of the corresponding control value to be considered as no significant interference.

No significant interference was observed when the following substances were tested for interference with this methodology. The data (calculated according to NGSP/DCCT) of interference studies on BS-400 are summarized below.

Interfering Substance	Interferent Concentration (mg/dL)	Analyte Concentration (HbA1c%)	Relative Deviation (%)*
Ascorbic acid	30	9.03	-2.21
Bilirubin	50	8.93	+2.61
Intralipid	2000	8.87	-2.26
Glucose	1000	8.97	-4.83

**Representative data, results in different instruments, laboratories may vary.*

■ Method Comparison

Correlation studies were performed using CLSI Approved Guideline EP09¹¹. The Mindray System (Mindray BS-800/Mindray HbA1c Reagent) (y) was compared with comparison system (Mindray H50/Mindray HbA1c Reagent) (x) using the same specimens. The statistical data (calculated according to NGSP/DCCT) obtained by linear regression are shown in the table below *:

Regression Fit	Correlation Coefficient (r)	Sample (N)	Concentration Range
$y=0.9645x+0.1424$	0.9967	183	3.8%-15.5%

**Representative data, results in different instruments, laboratories may vary.*

Result interpretation

The results could be affected by drugs, disease, or endogenous substances¹². When the reaction curve is abnormal, it is recommended to retest and check the result.

XX

HbA1c

These sections are applied for HbA1c assayed using whole blood samples without Centrifugation. The different part is as below.


Specimen collection and preparation

■ Specimen types

K_2 -EDTA, K_3 -EDTA, NaF-EDTA, lithium heparin and sodium heparin whole blood are suitable for samples.

■ Preparation for Analysis

- 1.Use the suitable tubes or collection containers and follow the instruction of the manufacturer; avoid effect of the materials of the tubes or other collection containers.
- 2.Do not centrifuge samples.
- 3.Refer to the flowchart below:

■ On-board lysis

1. Refer to "BS-series System Operations Manual" -> "General Operating Procedure" -> "Whole Blood Test".
2. Mix all specimens thoroughly by low speed vortexing or gently inverting 10 times prior to loading onto the Mindray BS series System.

- 3.The proportion of whole blood sample and Pretreatment Solution is 10ul:200ul.
- 4.Warning: For the anticoagulation tubes, the height of total sample in the tube should be into the range of 10mm-55mm. If out of the range, adjust to the height of sample or change to the manual lysis.
- 5.Warning: Micro-sample cups and Eppendorf microfuge tube are not allowed.

■ Manual lysis

- 1.Mix all specimens thoroughly by low speed vortexing or gently inverting 10 times.
- 2.Get 10 μ L whole blood into a sample cup or an Eppendorf microfuge tube. The whole blood is easily adsorbed to the sample probe and pipette, which may cause exceed the target volume.
- 3.Add 200 μ L of Pretreatment Solution into test tube.
- 4.Close the test tube and lysis the blood by inverting 10 times or shaking for 10s on the vortex mixer.
- 5.The hemolysate can be used as working samples after 4 minutes.
- 6.The volume of whole blood and Pretreatment Solution can be amplified in equal proportion if necessary.
- 7.Specimens should be tested as soon as possible after pre-analytical treatment.

■ Sample Stability

Whole blood Stability: 8 hours at 15-25°C

7 days at 2-8°C⁸

30 days at (-25)-(-15) °C

Hemolysate Stability: 8 hours at 15-25°C

24 hours at 2-8°C

Sample stability claims were established by manufacturer and/or based upon references, each laboratory should establish its own sample stability criteria.

Materials required but not provided

- 1.Mindray chemistry analyzers BS-600M or others which have the function of HbA1c assayed using whole blood samples without Centrifugation and General laboratory equipment.

Assay procedure

Twin test: The Hb and HbA1c concentration were test by using reagent R1/R2.

HbA1c

1.Hb

Parameters Item	BS-600M chemistry analyzers
Assay type	Endpoint
Wavelength (Primary/Secondary)	505/800 nm
Reaction direction	Increase
R1	90 μ L
incubate at 37°C for 1 min, read the absorbance A1, then add:	
Sample or Calibrator	8 μ L
Mix thoroughly, incubate at 37°C for 5 min, read the absorbance A2, Then calculate $\Delta A_{\text{Hb-calibrator}}=(A2-A1)$ calibrator, $\Delta A_{\text{Hb-sample}}=(A2-A1)$ sample	

2.HbA1c

Parameters Item	BS-600M chemistry analyzers
Assay type	Endpoint
Wavelength (Primary/Secondary)	660/800 nm
Reaction direction	Increase
R1	90 μ L
Sample or Calibrator	8 μ L
Mix, incubate at 37°C for 5 min, read the absorbance A1, then add:	
R2	30 μ L
Mix thoroughly, incubate at 37°C for 5 min, read the absorbance A2, Then calculate $\Delta A_{\text{HbA1c-calibrator}}=(A2-A1)$ calibrator, $\Delta A_{\text{HbA1c-sample}}=(A2 -A1)$ sample	

3.HbA1c%

HbA1c% is the most common form in the reports. Report results of hemoglobin A1c in mmol/mol (IFCC) or percent (NGSP):

Protocol 1 according to IFCC: $\text{HbA1c (mmol/mol)} = \text{HbA1c/Hb} \times 1000$

Protocol 2 according to NGSP/DCCT: $\text{HbA1c(\%)} = 91.5 \times \text{HbA1c /Hb} + 2.15$

Parameters may vary in different chemistry analyzers, may adjust in proportion if necessary. For Mindray BS series chemistry analyzers, Reagent Parameters is available on request. Please refer to the appropriate operation manual for the analyzers.

“Calibration” and “Quality control” refer to the section of “HbA1c assayed using whole blood samples with Centrifugation”.

Calculation

$C \text{ sample} = (\Delta A \text{ sample} / \Delta A \text{ calibrator}) \times C \text{ calibrator}$.

The BS series chemistry analyzer detects the change of absorbance(ΔA) and calculates the Hb/HbA1c concentration of each sample automatically after calibration.

The system calculates HbA1c and Hb results in $\mu\text{mol/L}$ (SI units). To convert the results to g/dL (common units), use the following equation:

Conversion factor: $\mu\text{mol/L} \times 0.1354 = \text{g/dL}$.

Dilution

If the value of sample exceeds $300 \mu\text{mol/L}$ (Hb) or $39 \mu\text{mol/L}$ (HbA1c), the sample should be diluted with Pretreatment Solution (e.g. 1+1) and rerun; the result of Hb and HbA1c should be multiplied by 2. No conversion factor is required for the mmol/mol HbA1c (IFCC) and % HbA1c (DCCT/NGSP) result.

Expected values⁹

	Sample Type	S.I. Units
Whole blood	According to IFCC	20.20-42.06 mmol/mol
	According to NGSP/DCCT	4.0%-6.0%

The expected value is provided from reference, and Mindray has verified it by 60 samples of people from southern China.

Each laboratory should establish its own reference intervals based upon its particular locale and population characteristics since expected values may vary with geography, race, sex and age.

Performance Characteristics

■ Analytical Sensitivity

The Limit of Blank (LoB) and Limit of detection (LoD) were determined on BS-600M in accordance with CLSI Approved Guideline EP17-A2.

Hb:	HbA1c:	HbA1c%(NGSP/DCCT):
LoB=4.55 $\mu\text{mol/L}$	LoB=1.12 $\mu\text{mol/L}$	LoB=2.97%
LoD=4.94 $\mu\text{mol/L}$	LoD=1.17 $\mu\text{mol/L}$	LoD=3.01%

The Limit of Blank is the 95th percentile value from $n \geq 60$ measurements of one or several analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95%.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples. The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of Blank with a probability of 95%).

■ Measuring range

The Mindray BS series systems provide the following linearity range:

Sample Type	S.I. Units	
	HbA1c%	3%-16% (NGSP/DCCT)
Whole blood	HbA1c	2-39 $\mu\text{mol/L}$
	Hb	60-300 $\mu\text{mol/L}$

HbA1c%: A high concentration HbA1c% sample (approximately 16%) is mixed with a low concentration sample (<3%) at different ratios, generating a series of dilutions. The HbA1c% concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient $r \geq 0.9900$.

HbA1c (Hb concentration $\geq 100 \mu\text{mol/L}$): A high concentration HbA1c sample (approximately 39 $\mu\text{mol/L}$) is mixed with a low concentration sample (<2 $\mu\text{mol/L}$) at different ratios, generating a series of dilutions. The HbA1c concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient $r \geq 0.9900$.

Hb: A high concentration Hb sample (approximately 300 $\mu\text{mol/L}$) is mixed with a low concentration sample (<60 $\mu\text{mol/L}$) at different ratios, generating a series of dilutions. The Hb concentration of each dilution is determined using Mindray System, the linearity range is demonstrated with the correlation coefficient $r \geq 0.9900$.

■ Precision

Precision was determined by following CLSI Approved Guideline EP05-A3¹⁰, each sample was assayed 2 times per run, 2 runs per day, for a total 20 days. The precision data (calculated according to NGSP/DCCT) of controls and human samples on BS-600M are summarized below*.

Specimen Type (N=80)	Mean (%)	Repeatability		Within-Lab	
		SD (%)	CV %	SD (%)	CV %
Control N	5.16	0.02	0.35	0.02	0.44
Contro P	10.28	0.02	0.18	0.03	0.29
Whole blood 1	5.47	0.02	0.32	0.03	0.54
Whole blood 2	6.63	0.02	0.33	0.03	0.45
Whole blood 3	9.32	0.02	0.22	0.04	0.46
Whole blood 4	11.38	0.03	0.27	0.05	0.41

*Representative, results in different instruments, laboratories may vary.

■ Analytical Specificity

1. Endogenous substances and Durgs

The samples with different concentration interfering substance were prepared by addition of interferent to human sample pools, and recovers are within $\pm 5.0\%$ of the corresponding control value to be considered as no significant interference.

No significant interference was observed when the following substances were tested for interference with this methodology. The data (calculated according to NGSP/DCCT) of interference studies on BS-600M are summarized below*.

Interfering Substance	Interferent Concentration	Analyte Concentration	Relative Deviation
Ascorbic acid	3 mg/dL	6.64%	+0.41%
Conjugated bilirubin	15 mg/dL	6.63%	-2.06%
Unconjugated bilirubin	10 mg/dL	6.66%	-2.19%
Triglycerides	2000 mg/dL	6.67%	-0.10%
Melbine	5.1 mg/dL	6.61%	-0.40%
Ibuprofen	50 mg/dL	6.59%	+0.15%
Acarbose	50 mg/dL	6.58%	-0.08%
Panadol	200 μ g/mL	6.59%	-3.07%
Aspirin	50 mg/dL	6.56%	0.00%

*Representative data, results in different instruments, laboratories may vary.

2. Hemoglobin Derivatives

The samples with different type of hemoglobin derivatives were obtained by addition of potential interferent to human sample pools, and recovers are within $\pm 5.0\%$ of the corresponding control value to be considered as no significant interference.

No significant interference was observed when the following substances were added with this methodology. The data (calculated according to NGSP/DCCT) of interference studies on BS-600M are summarized below*.

Type of Hemoglobin Derivatives	Analyte Concentration	Relative Deviation
Acetylated Hemoglobin generated by adding ≥ 50 mg/dL of Aspirin	6.50%	+0.71%
Carbamylated Hemoglobin generated by adding ≥ 10 mmol/L of Cyanate	6.61%	-1.43%
Labile Hemoglobin generated by adding ≥ 1000 mg/dL of Glucose	6.61%	-1.13%

*Representative data, results in different instruments, laboratories may vary.

3. Hemoglobin Variants

Specificity was assessed by comparing the Hemoglobin A1c values to reference values for samples containing abnormal hemoglobins, and recovers are within $\pm 5.0\%$ of the corresponding control value to be considered as no significant interference.

Heterozygous hemoglobin variants (Hb-C, Hb-D, Hb-E, Hb-S, HbA2) do not interfere with the Hemoglobin A1c assay.

The statistical data (calculated according to NGSP/DCCT) of hemoglobin variants study on BS-600M are summarized in the following table.

Hemoglobin Variant Relative % Difference from Reference Concentration

Variant	~ 6.5 %HbA1c	~ 9.0 %HbA1c
Hb-C	-3.42%	-3.29%
Hb-D	-1.24%	-0.74%
Hb-E	+1.53%	+1.36%
Hb-S	-1.05%	-1.39%
HbA2	+3.22%	+3.85%
Hb-F	Difference exceeds -5% when the amount of HbF in the sample exceeds 5%*.	

*A negative % difference with HbF is proportional in magnitude to the % HbF present in the sample. For example, when the amount of HbF in the sample was 17%, the % difference was -11.6% on the Mindray BS-600M System.

NOTE: The presence of multiple variants in a sample may impact the % difference.

■ Method Comparison

Correlation studies were performed using CLSI Approved Guideline EP09-A3¹¹. The Mindray System (Mindray BS-2800M/Mindray HbA1c Reagent) (y) was compared with comparison system (Mindray H50/ Mindray HbA1c Reagent) (x) using the same specimens. The statistical data (calculated according to NGSP/DCCT) obtained by linear regression are shown in the table below *:

Regression Fit	Correlation Coefficient (r)	Sample (N)	Concentration Range
$y = 1.0178x - 0.08688$	0.9962	118	4.72%-15.65%

*Representative data, results in different instruments, laboratories may vary.

Result Interpretation

The results could be affected by drugs, disease, or endogenous substances¹². When the reaction curve is abnormal, it is recommended to retest and check the result.

XX

Warnings and precaution

1. For in vitro diagnostic use only. For laboratory professional use.
2. The test is designed only for accurate and precise measurement of HbA1c%. The individual results for total Hb and HbA1c concentration should not be reported.
3. Please take the necessary precautions for handling all laboratory reagents.
4. Please confirm the integrity of the package before use. Do not use the kits with damaged packages. The reagents avoid direct exposure to sunlight and freezing. The results can't be assured when stored at inappropriate condition.
5. If unintentionally opened before used, store the reagents tightly capped at 2-8°C and protected from light, and the stability is equally to in-use stability.
6. Do not mix reagents with different lots.
Do not use the reagents beyond the expiration date and the in-use date.
Do not mix fresh reagents with in-use reagents.
Avoid the formation of foam.
7. Instability or deterioration should be suspected if there are visible signs of leakage, precipitates or microbial growth, or if calibration or controls do not meet the insert and/or the Mindray System criteria.
8. Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
9. Preservative contained. Do not swallow. Avoid contact with skin and mucous membranes.
10. When the reagents accidentally enter the eyes and mouth, or contact with the skin, immediately wash with plenty of water. If necessary, visit the doctor for further medical treatment.
11. Safety data sheet is available for professional user on request.
12. Disposal of all waste material should be in accordance with local guidelines.
13. All human material should be considered potentially infectious.
14. All identified risks have been reduced as far as possible without adversely affecting the benefit-risk ratio, and the overall residual risk is acceptable.
15. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
16. This kit contains components classified as follows in accordance with the Regulation (EC) No 1272/2008:

**Danger**

H317	May cause an allergic skin reaction.
H360D	May damage the unborn child.
H412	Harmful to aquatic life with long lasting effects.

Prevention:

P201	Obtain special instructions before use.
P280	Wear protective gloves and protective clothing.
P261	Avoid breathing mist/vapours/spray.
P273	Avoid release to the environment.
P272	Contaminated work clothing should not be allowed out of the workplace.

Response:

P308+P313	IF exposed or concerned: Get medical advice/attention.
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before reuse.

Disposal:

P501	Dispose of contents/container to authorised hazardous or special waste collection point in accordance with any local regulation.
------	--

References

1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;1441-1447.
2. Goldstein DE, Little RR, Lorenz RA et al. Tests of glycemia in diabetes. Diabetes Care, 1995;18:896-909.
3. Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. Diabetes Care, 1994;17:938-939.
4. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med, 1993;329:977-986.

5. Santiago JV. Lessons from the diabetes control and complications trial. *Diabetes*, 1993;42:1549-1554.
6. UK Prospective Diabetes Study (UKPDS) group. Intensive blood-glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). *Lancet*, 1998;352:837-853.
7. Nathan DM, Kuenen J, Borg R et al. Translating the A1c Assay Into Estimated Average Glucose Values. *Diabetes Care*, 2008; 31:1473-1478.
8. Wu, Alan HB. *Tietz clinical guide to laboratory tests*. 4th ed. Elsevier Health Sciences, 2006; 480.
9. Hanas R, John G, International HbA1c Consensus Committee. 2010 Consensus statement on the worldwide standardization of the hemoglobin A1C measurement. *ClinChem Lab Med*, 2010;48:775-776.
10. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
11. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition*. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
12. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington, DC: AACC Press; 2000:3-433,3-434.

Graphical symbols

In Vitro Diagnostic
medical device

Unique device
identifier

European
Conformity

Consult Instructions
For use

Use-by
date

Authorized representative in
the European Community

Batch Code

Temperature
limit

Manufacturer

Catalogue
number

Keep away from sunlight

Indicates a medical device that needs protection from light sources

HbA1c

mindray

© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998; **Fax:** +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;1441-1447.
2. Goldstein DE, Little RR, Lorenz RA et al. Tests of glycemia in diabetes. Diabetes Care, 1995;18:896-909.
3. Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. Diabetes Care, 1994;17:938-939.
4. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med, 1993;329:977-986.
5. Santiago JV. Lessons from the diabetes control and complications trial. Diabetes, 1993;42:1549-1554.
6. UK Prospective Diabetes Study (UKPDS) group. Intensive blood-glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). Lancet, 1998;352:837-853.
7. Nathan DM, Kuenen J, Borg R et al. Translating the A1c Assay Into Estimated Average Glucose Values. Diabetes Care, 2008; 31:1473-1478.
8. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 480.

- Hanas R, John G, International HbA1c Consensus Committee. 2010 Consensus statement on the worldwide standardization of the hemoglobin A1C measurement. ClinChem Lab Med, 2010;48:775-776.
- CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
- Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-433,3-434.

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Адрес: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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Referências

1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;1441-1447.
2. Goldstein DE, Little RR, Lorenz RA et al. Tests of glycemia in diabetes. *Diabetes Care*, 1995;18:896-909.
3. Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. *Diabetes Care*, 1994;17:938-939.
4. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med*, 1993;329:977-986.
5. Santiago JV. Lessons from the diabetes control and complications trial. *Diabetes*, 1993;42:1549-1554.
6. UK Prospective Diabetes Study (UKPDS) group. Intensive blood-glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). *Lancet*, 1998;352:837-853.
7. Nathan DM, Kuenen J, Borg R et al. Translating the A1c Assay Into Estimated Average Glucose Values. *Diabetes Care*, 2008; 31:1473-1478.
8. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 480.
9. Hanas R, John G, International HbA1c Consensus Committee. 2010 Consensus statement on the worldwide standardization of the hemoglobin A1C measurement. *ClinChem Lab Med*, 2010;48:775-776.
10. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
11. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
12. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-433,3-434.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Morada: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. China

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3. Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. *Diabetes Care*, 1994; 17:938-939.
4. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med*, 1993; 329:977-986.
5. Santiago JV. Lessons from the diabetes control and complications trial. *Diabetes*, 1993; 42:1549-1554.
6. UK Prospective Diabetes Study (UKPDS) group. Intensive blood-glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). *Lancet*, 1998; 352:837-853.
7. Nathan DM, Kuenen J, Borg R et al. Translating the A1c Assay Into Estimated Average Glucose Values. *Diabetes Care*, 2008; 31:1473-1478.
8. Wu, Alan HB. *Tietz clinical guide to laboratory tests*. 4.^a edición Elsevier Health Sciences, 2006; 480.
9. Hanas R, John G, International HbA1c Consensus Committee. 2010 Consensus statement on the worldwide standardization of the hemoglobin A1C measurement. *ClinChem Lab Med*, 2010; 48:775-776.
10. CLSI. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition*. documento EP05-A3 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
11. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition*. Documento EP09-A3 del CLSI. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
12. Young DS. *Effects of Drugs on Clinical Laboratory Tests*, 5.^a edición Washington, DC: AACC Press; 2000:3-433, 3-434.

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Bibliografia

1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;1441-1447.
2. Goldstein DE, Little RR, Lorenz RA et al. Tests of glycemia in diabetes. Diabetes Care, 1995;18:896-909.
3. Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. Diabetes Care, 1994;17:938-939.
4. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med, 1993;329:977-986.
5. Santiago JV. Lessons from the diabetes control and complications trial. Diabetes, 1993;42:1549-1554.
6. UK Prospective Diabetes Study (UKPDS) group. Intensive blood-glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). Lancet, 1998;352:837-853.
7. Nathan DM, Kuenen J, Borg R et al. Translating the A1c Assay Into Estimated Average Glucose Values. Diabetes Care, 2008; 31:1473-1478.
8. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 480.
9. Hanas R, John G, International HbA1c Consensus Committee. 2010 Consensus statement on the worldwide standardization of the hemoglobin A1C measurement. ClinChem Lab Med, 2010;48:775-776.
10. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
11. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
12. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-433,3-434.

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Indirizzo: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, Repubblica popolare cinese

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Referanslar

1. Burtis CA, Ashwood ER, Bruns DE, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 5th ed. Elsevier Saunders 2012;1441-1447.
2. Goldstein DE, Little RR, Lorenz RA et al. Tests of glycemia in diabetes. Diabetes Care, 1995;18:896-909.
3. Goldstein DE, Little RR. More than you ever wanted to know (but need to know) about glycohemoglobin testing. Diabetes Care, 1994;17:938-939.
4. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med, 1993;329:977-986.
5. Santiago JV. Lessons from the diabetes control and complications trial. Diabetes, 1993;42:1549-1554.
6. UK Prospective Diabetes Study (UKPDS) group. Intensive blood-glucose control with sulfonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). Lancet, 1998;352:837-853.
7. Nathan DM, Kuenen J, Borg R et al. Translating the A1c Assay Into Estimated Average Glucose Values. Diabetes Care, 2008; 31:1473-1478.
8. Wu, Alan HB. Tietz clinical guide to laboratory tests. 4th ed. Elsevier Health Sciences, 2006; 480.
9. Hanas R, John G, International HbA1c Consensus Committee. 2010 Consensus statement on the worldwide standardization of the hemoglobin A1C measurement. ClinChem Lab Med, 2010;48:775-776.
10. CLSI. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

11. CLSI. Measurement Procedure Comparison and Bias Estimation Using Patient Sample; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
12. Young DS. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington, DC: AACC Press; 2000:3-433,3-434.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, High-Tech Industrial Park, Nanshan, Shenzhen, 518057, P. R. Çin

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© 2022-2023 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés.
Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)