INgezim® West Nile COMPAC

Prod Ref: 10.WNV.K3

Ensayo inmunoenzimático de bloqueo para la detección de anticuerpos específicos del West Nile Virus.
Blocking immunoenzymatic assay for the specific detection of antibodies to West Nile Virus.

Version Doc:27.05.2024

Nº de registro en España: 1815-RD

Registration number in Spain: 1815-RD

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KIT COMPOSITION

<i>REAGENT</i>	<i>2 Plates</i>		<i>5 Plates</i>	
	Uni.	Vol.	Uni.	Vol.
Microtitration plates	2		5	
Vials with Positive Control, ready to use	1	1,25mL	2	1,25mL
Vials with Negative Control, ready to use	1	1,25mL	2	1,25mL
Vials with Peroxidase Conjugate, ready to use	1	15mL	2	15mL
Bottles with Washing Solution 25x concentrated	1	125mL	1	125mL
Bottles with Diluent (DE01-01), ready to use	1	125mL	1	125mL
Bottles with Substrate (TMB), ready to use	1	30mL	1	60mL
Bottles with Stop Solution, ready to use	1	60mL	1	60mL

OTHER MATERIALS AND REAGENTS NEEDED NOT PROVIDED WITH THE KIT

Distilled or deionized water.
Micropipettes from 5 to 200 µL.
Disposable micropipette tips.
Washing plates device.
Test tubes from 50 to 250 mL.
ELISA Reader (450 nm filter)

TECHNICAL BASIS

The kit has been designed to detect antibodies specific for West Nile virus. The kit is able to detect a very low titer of antibodies in sera of different infected animal species (birds, horses, humans, etc.).

The INgezim® West Nile COMPAC kit is based on the blocking enzyme immunoassay described below:

The plates are coated with inactivated viral antigen. After adding the sample to the well, if it contains specific antibodies against WNV, they will bind to the antigen absorbed on plate while if the sample does not contain specific antibodies, they will not bind to the antigen.

If we add a specific monoclonal antibody (conjugated with peroxidase) against the ectodomain of glycoprotein E of WNV, it will compete with the antibodies of the serum. If the serum samples contain specific antibodies, they will not permit the binding of the labelled Mab to the antigen whereas if it does not contain specific antibodies the Mab will bind to the antigen on the plate. After washing the plate to eliminate all non-fixed material, we can detect the presence or absence of labelled Mab by adding the substrate (TMB) that in presence of the peroxidase will develop a colorimetric reaction.

PRECAUTIONS AND WARNINGS FOR USERS

1. Read the instructions of use carefully.
2. Bring all reagents to room temperature prior to use.
3. Do not mix instructions or reagents from different kits and avoid any contamination of the reagents.
4. Do not use components after expiration dates and do not mix components from different lots.
5. There should be no eating, drinking, or smoking where specimens or Kit reagents are being handled and do not pipette by mouth.
6. Use a new tip for each serum sample.
7. The substrate is extremely sensitive to light and contamination, so it is recommended to remove only the volume to be used and never return the excess substrate to the bottle.
8. Stop solution is a strong acid solution. In case of accidental contact with skin, wash gently with water.

STORAGE OF COMPONENTS

Keep all the reagents between +2°C and +8°C.

INFORMATION ABOUT THE WASHING STEPS

The washing steps could be done using an automatic washing machine or a multichannel pipetting device suitable for dispensing 300 µL on each well. After the incubation periods, the washing steps must be done following the next instructions:

- Throw out the content of the plate by a brusque turnover of the plate to avoid the possible mixture of the content from one well to another. As a precaution, the wells should be emptied over a cuvette containing 0.1M NaOH solution, as the test samples may contain infectious agents.
- Dispense a volume of 300 µL of washing solution on each well.
- Shake delicately the plate, avoiding the contamination between wells.
- Turn over the plate brusquely to empty the wells.
- Repeat the process as much times as is indicated on the instructions of the Kit. After the washing step, verify that the next reagent to be added to the plate is ready to use. Do not maintain the plate on dry.
- After the last step of washing shake the plate turned over an absorbent filter paper.

PREPARATION OF SAMPLES

Sera samples must be tested at 1/5 dilution. This dilution could be made directly on the assay plate by adding 40 µL of diluent and 10 µL of sample to each well.

Mix gently to homogenate the dilution in the well.

PREPARATION OF REAGENTS

Washing solution:

Dilute one part of the concentrate washing solution provided in the Kit with 24 parts of distilled or deionized water (40 mL of the concentrated solution in 960 mL of distilled or deionized water). Once ready this solution remains stable between +2°C and +8°C.

Control sera:

They are supplied ready to use. **Do not dilute.**

Preparation of the conjugate:

Conjugate is ready to use and must be used adding 50 µL/well. **Do not dilute.**

TEST PROCEDURE

1. All reagents must be allowed to come to room temperature before use.
2. Add samples as indicated in section "Preparation of Samples". Add 50 µL of controls in duplicate wells. Gently mix the content of the wells, seal the plate and incubate overnight (**16-24 hours**) at **room temperature (20°C to 25°C)**.
3. Wash 3 times following the described procedure.
4. Add 50 µL of conjugate ready to use to each well. Seal the plate and **incubate for 1h at room temperature.**
5. Wash 5 times following the described procedure.
6. Add 50 µL of substrate, to each well. Keep the plate for **15 min at room temperature.**
7. Add 50 µL of stop solution to each well.
8. Read the OD of each well at 450 nm within 5 min after the addition of stop solution.

READING AND RESULT INTERPRETATION

In case that the samples had been run in duplicate, it has to be considered the mean of both OD values. In the same way, you have to make the mean of the values obtained in the two wells of positive control and the two wells of negative control.

VALIDATION OF THE RESULTS:

- The OD in the wells with Negative Control (NC) must be > 0.8
- The OD in the wells with Positive Control (PC) must be < 0.35

INTERPRETATION OF THE RESULTS:

Calculate the inhibition percentage (IP) of each sample as follow:

$$\text{IP of sample} = 100 - \left(\left[\frac{\text{Sample OD}}{\text{NC OD}} \right] \times 100 \right)$$

- Samples will be considered **positive** (there are antibodies specific for WNV), when the IP is ≥ 40%
- Samples will be considered **negative** (there are no antibodies specific for WNV in the sample) when the IP is ≤ 30%
- Samples with IP between 30% and 40% should be considered **doubtful** and retested by seroneutralization (SN) assay.

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Designed & manufactured in Spain by:

GOLD STANDARD DIAGNOSTICS MADRID SA
C/Hermanos García Noblejas, 41 2ª planta

28037- MADRID (SPAIN)

Tlf: +34 91368.05.01/04

Fax: +34 91 408.75.98

E-mail: info.spain@eu.goldstandarddiagnostics.com

www.goldstandarddiagnostics.com

Distribuido en
Distributed in

por
by