

Digital Angiography System

# DAR-9500f







## Operation Manual

Read the instruction manual thoroughly before you use the product.  
Keep this instruction manual for future reference.

---

## Categories of Information

In this manual, safety and utilization information is categorized and indicated as follows.

| Symbol   | Significance  |
|--|---|
|  <b>DANGER</b>  | Indicates critically hazardous situation that, if not avoided, may result in serious injury or death.                       |
|  <b>WARNING</b> | Indicates indirectly or potentially hazardous situation that, if not avoided, may result in serious injury or death.        |
|  <b>CAUTION</b> | Indicates hazardous situation that, if not avoided, may result in minor or moderate injury, damage to the product, or fire. |
|  <b>NOTE</b>   | Indicates information for proper use of the product.  |
|               | Indicates information for convenient use of the product.  |
|               | Indicates reference for additional information.   |

# Introduction

This operator's manual describes operation of the DAR-9500f Digital Angiography System.

Please read this manual thoroughly before using this equipment.

If precautions described in this manual are not kept, it may cause damage to the system, operators and patients. And also, it is difficult to describe an unpredictable precautions, so please contact our service representatives if you use the system by not described procedure or have any questions in this manual.

This manual should be kept available for future reference.

////////////////////////////////////

## Disclaimer

SHIMADZU Corporation reserves the copyright of this Manual and no part of its content may be copied or reproduced without the company's permission.

The contents of this Manual may be changed for improvement without notice.

Despite commitment and effort, errors and omissions found later in this Manual may not be immediately corrected in some cases.

Windows is either a registered trademark or trademark of Microsoft Corporation in the United States and/or other countries.

SHIMADZU Corporation is not liable for any consequence of operation following this Manual.

////////////////////////////////////

---

## Precautions for Installation and Use of This System

Please observe the following “Operating Precautions for Safety in the Use of Electric Medical Equipment”.

### **Operating Precautions for Safety in the Use of Electric Medical Equipment**

1. Only an experienced technician should operate the equipment.
2. When installing the equipment, pay attention to the following items:
  - (1) Do not install system near water faucet or similar equipment.
  - (2) Install it away from potential sources of problems such as abnormal pressure, temperature or humidity, drafts, direct sunlight, chlorine dust or sulfur gas.
  - (3) During transportation and operation of the equipment, avoid tilting, vibration and any impact.
  - (4) Keep the equipment away from areas where chemicals or gases are stored.
  - (5) Use only the correct electrical power source with matching frequency, voltage and current (or wattage).
  - (6) Check the condition of the battery power source (power and polarity) before operating the equipment.
  - (7) Properly ground the equipment.
3. Before operating the equipment, pay attention to the following items:
  - (1) Check the conditions of switch contacts, polarity, dial settings, and meters, and make sure the equipment performs correctly.
  - (2) Confirm that the ground is connected properly.
  - (3) Check all wiring for proper and correct connections.
  - (4) Pay attention when using more than one unit at a time, because it may lead to an incorrect diagnosis and cause complications.
  - (5) Check the conditions of any external electric circuit, that will be directly connected to a patient.
  - (6) Check the condition of the battery power source.



- 
4. While operating the equipment, pay attention to the following items:
    - (1) Do not exceed recommended times or the amount of radiation needed for diagnosis or therapy.
    - (2) Observe the equipment and patient continuously for early detection of problems.
    - (3) When a problem is detected with the equipment, take proper action to stop the equipment without harming the patient.
    - (4) Do not let the equipment touch the patient.
  5. After operating the equipment, pay attention to the following items:
    - (1) Turn off the switches and return the dials to their original position in the prescribed order. Then, turn off the main power switch.
    - (2) Do not pull the power cable forcibly from the outlet using the power cord.
    - (3) When storing the equipment, pay attention to the following factors:
      - (i) Keep it away from the water.
      - (ii) Store it away from the potential causes of problems such as abnormal pressure, temperature or humidity, draft, direct sunlight, chlorine dust or sulfur gas.
      - (iii) Avoid tilting, vibration and any impact when storing.
      - (iv) Store the equipment away from areas where chemicals and gases are stored.
    - (4) Clean all attachments, cables and contacts, and store them in one place.
    - (5) Keep the equipment clean to avoid problems during the next use.
  6. When the equipment is found to be out of order, do not try to repair it. Immediately call a certified repair technician for repair.
  7. Do not modify any part of the equipment.
  8. Preventive maintenance:
    - (1) The equipment and its parts should be periodically checked.
    - (2) If the equipment has not been in operation for an extended period of time, test it prior to actual operation to make sure it works correctly and safely before use.
  9. Operate properly according to the operating manual.

## Precautions in Use



### WARNING

**The responsibility for managing use and maintenance of medical equipment lies with the user.**

Use of this device is restricted to a diagnostic radiology technician or a person with a certificate indicating equal proficiency.



### WARNING

**Do not modify the equipment.**

Repair and inspection of the inside of the equipment is dangerous. Be sure to contact our service agency for repair and inspection of any kind.



### WARNING

**Perform periodical maintenance.**

Maintenance is required to assure the safety and performance of this equipment. For details on the maintenance inspections to periodically be performed by the operator, please refer to the descriptions contained in this operation manual.



### WARNING

**Repair and maintenance of this equipment can only be performed by engineers approved by SHIMADZU.**

Maintenance must be assigned to specially trained experts. Contact the Shimadzu Service Representative for repair and maintenance.

**WARNING****Prepare alternative system.**

Prepare alternative system with fluoroscopy function, in case any trouble happens and fluoroscopy/radiography can not be operated normally.

**WARNING****Beware of X-ray exposure.**

Improper use of the X-ray equipment might cause the operator or patient to be accidentally exposed to X-ray radiation. During X-ray radiation, any person other than the subject patient should not stay in the irradiation room. If circumstances compel any non-subject person to stay in the room, ample protection should be provided for that person.

**WARNING****Do not splash water on the equipment.**

Splashing water might cause an electric shock. When cleaning the equipment, use a cloth dipped in an antiseptic solution (Medical Alcohol) and squeeze out all excess liquid before wiping only the surface.

**WARNING****Accuracy of displayed values are not guaranteed.**

Displayed values measured by the measurement functions of this equipment are not absolute values but relative values based on the capability of the instruments used.

**WARNING**

**Do not use the Injector which is not conformed to the standard of IEC60601/1.**



## WARNING

**If the operator has no experience in operating the equipment, be sure that he or she receives instruction on how to operate it from our engineers or someone who has enough experience to use the equipment.**

In order to operate the equipment, safely, an explanation of the operation needs to be lectured. When installing the equipment, our engineers explain the operating procedure. Follow their directions and operate the equipment correctly.



## WARNING

**Do NOT perform any maintenance work on any part of the equipment during clinical use.**

It may cause injury.



## CAUTION

For California, USA Only

This product contains a battery that contains perchlorate material.

Perchlorate Material-special handling may apply.


See [www.dtsc.ca.gov/hazardouswaste/perchlorate](http://www.dtsc.ca.gov/hazardouswaste/perchlorate)

## Connector



### WARNING

Do not connect other devices which are not specified by SHIMADZU, to any of the connectors (enclosed in a dashed-line box in the illustration below).



## Other Software Installation



### WARNING

Do not install other software on the computer. Do not connect peripheral equipment, or modify the computer hardware. Otherwise, the system may fail to start because of a change in the OS environment or a driver conflict. Not heeding this warning will invalidate the warranty.

## Internet Connection



### CAUTION

Do not connect this system to the Internet. Shimadzu shall not guarantee the security against virus infection or leak of hospital/patient information through Internet connection. Do not connect this system to other computer networks that can be linked to the Internet.

### ■ Scan Converter BNC Cable Positioning



#### CAUTION

*(Applicable only to the BNC cables provided with the Scan Converter option.)* Do not run the scan converter BNC video cables next to power or signal cables because this may cause video noise to appear.

### ■ Data Erasure or Loss



#### CAUTION

Data saved on magnetic disk may be erased because of improper operation or system accident. Be sure to backup important data to external storage such as CD-R.

Since it is possible that external backups may become damaged, multiple backups are recommended.

Note that Shimadzu will not be held responsible for data erasure or loss.

### ■ Disposal

If you must dispose of the product, the following caution must be observed:



#### CAUTION

The device contains substances which may pollute the environment if disposed of incorrectly. Contact our service office or representative for disposal of the system or any parts.

**Cautions on Wireless Communication Using Mouse****WARNING**

This system operates in the same frequency band as premises radio stations for identifying mobile devices in factory production lines (license required) and specified low-power radio stations (license not required) in addition to industrial devices (ex. microwave ovens), scientific instruments, and medical equipment complies with CISPR EMISSION requirements.

The use of this system may result in RF interference with the above mentioned equipment and radio stations. Make sure you understand and heed the following cautions when operating the system.

- Before using this system, make sure that no RFID premises radio stations and specified low-power radio stations or similar equipment are used in the immediate area.
- If this system causes damaging RF interference to affect an RFID premises radio station, stop using the system immediately and contact your Shimadzu service representative.
- If this system causes damaging RF interference to affect an RFID specified low-power radio station or amateur radio station, contact your Shimadzu service representative.
- This system may suffer interference from other equipment that emits radio waves (such as microwave ovens, Bluetooth devices, and digital cordless telephones). Use the system after moving such devices as far away as possible to prevent interference.

**CAUTION**

**Do NOT use frequency channel (2.4 GHz band) of wireless mouse for other wireless devices.**

## Network Connection



### CAUTION

1. Connect to the network in order to achieve the following purpose.
  - Acquire patient information, study information and past study image from the external equipment.
  - Send study image, study record, dose information and system maintenance information to the external equipment.
2. Connect to the IPv4 network that supports 1000 Base-T/100 Base-T.
3. Use UTP cable with enhanced category 5 (Cat 5e) or higher for LAN cable.
4. Intended information flow

| Destination               | Typical Information  |
|---------------------------|--|
| Printer                   | Study image data   |
| MWM Server<br>PPS Manager | Patient information, Study information<br>Study record, Dose information |
| DICOM Image Server        | Study image data   |
| Remote Maintenance Server | System maintenance information   |

5. Hazardous conditions in case of network failure


| Destination               | Hazardous Situation   |
|---------------------------|---|
| Printer                   | Cannot print study image.   |
| MWM Server<br>PPS Manager | Cannot acquire patient information and study information.<br>Cannot send study record and dose information. |
| DICOM Image Server        | Cannot send study image.<br>Cannot acquire past study image.  |
| Remote Maintenance Server | Cannot send system maintenance information.   |



## ■ Action for Environment

To all users of Shimadzu equipment in the European Union:

Equipment marked with the following symbol indicates that it was sold on or after 13th August 2005, which means it should not be disposed of with general household waste. Note that our equipment is for industrial/professional use only.



WEEE Mark

Contact a Shimadzu service representative when the equipment has reached the end of its life. They will advise you regarding the equipment take-back.

With your cooperation, we are aiming to reduce contamination from waste electronic and electrical equipment and preserve natural resources through re-use and recycling.

Do not hesitate to ask a Shimadzu service representative, if you require further information.

# Pediatric Imaging

## ■ General

X-ray diagnosis is convenient and useful to visualize inside of human body, however the use of X-ray has potential risk. This section describes pediatric imaging consistent with ALARA<sup>1</sup> principles to minimize the radiation dose.

Exposure to ionizing radiation is of particular concern in pediatric patients for three reasons:<sup>2</sup>

1. 1) younger patients are more radiosensitive than adults (i.e., the cancer risk per unit dose of ionizing radiation is higher for younger patients);
2. younger patients have a longer expected lifetime for the effects of radiation exposure to manifest as cancer; and
3. use of equipment and exposure settings designed for adult use can result in excessive radiation exposure for the smaller patient.

Besides this section, the some websites provide additional information for pediatric imaging.

American College of Radiology: <http://www.acr.org>

Image Gently: <http://www.imagegently.org>

US Food and Drug Administration, FDA: <http://www.fda.gov>

<sup>1</sup> ALARA: As Low As Reasonably Achievable

<sup>2</sup> FDA DRAFT GUIDANCE, Pediatric Information for X-ray Imaging Device Premarket Notifications, May 10, 2012

## ■ Suggestions for Pediatric Imaging

This section describes some suggestions for reducing dose of pediatric imaging. They are:

- minimizing unnecessary dose;
- adjusting parameters and
- monitoring dose indications.

### Minimizing unnecessary dose

To minimize unnecessary dose, adjust the collimator to cover just region concerned and the region concerned to be center of the image. Before exposure, use proper protective means in accordance with predefined guidelines, ex. protective clothes, when applicable, and check the parameters or the selection of proper preset which influence the radiation dose, based on ALARA principles.

### Adjusting parameters

Select a preset based on patient size or region concerned. That is a start point for adjusting parameters like kVp and MAS. Some of the systems equip detachable grid to reduce dose. Consultation with professionals to adjust proper parameters based on ALARA principles is strongly recommended.

### **Monitoring dose indications**

It is important to measure and quantify radiation dose for each patient. Some of the X-ray Diagnostic Imaging Systems equip the means to display estimates of entrance dose (dosearea product for general radiographic system or air kerma for interventional or noninterventional fluoroscopic system). Some also have means to store radiation dose information in the DICOM header of each image or means to generate separate radiation dose information besides clinical images. It is suggested to utilize those capabilities to monitor radiation dose, and consequently to minimize radiation dose to the patient.

## The Limited Product Warranty

The system warranty is for one year from the date of purchase.

The following failures or damage are not covered by the warranty:

1. Failure or damage due to the installation, relocation, or service not performed by a SHIMADZU Service Representative or a SHIMADZU designated contractor.
2. Failure or damage caused by products from other companies (except those purchased from SHIMADZU).
3. Failure or damage due to repairs using non-SHIMADZU certified service parts.
4. Failure or damage caused from not following the notices and procedures described in this manual.
5. Failure or damage due to an operating environment that is outside the requirements stipulated in this manual.
6. Failure or damage due to natural disasters such as fire, earthquake, flood, and lightning.

Service after the expiration of the warranty is available at a reasonable cost and should be performed by the SHIMADZU Service Representative.

## Revision History

| REVISION | DATE    | COMMENT   |
|----------|---------|---|
| Original | 2011.3  |   |
| A        | 2011.5  |   |
| B        | 2012.11 | <ul style="list-style-type: none"> <li>• Update cable and accessory lists.</li> <li>• Add SuperCine.</li> <li>• Add Pixel Spacing.</li> <li>• Add E-Shutters.</li> <li>• Delete Injection Delay Settings and add Acquisition Parameter Settings.</li> <li>• Add status of an image. (LIVE/LIH/STORED)</li> <li>• Add SCORE 3D (DA/DSA/CT).</li> </ul>   |
| D        | 2013.3  | <ul style="list-style-type: none"> <li>• Apply to Bi-plane system.</li> <li>• Add description about license information and dose report.</li> <li>• Update window images.</li> <li>• Update the Direct X.</li> <li>• Add adjustment procedure of brightness and contrast with mouse.</li> <li>• Add correction of image comment.</li> <li>• Update the name of each button and item.</li> </ul> |
| E        | 2013.7  | <ul style="list-style-type: none"> <li>• Add Peak Hold to DSA tool.</li> <li>• Update various images.</li> <li>• Update various buttons and items.</li> <li>• Add a NOTE on calibration after the power failure.</li> <li>• Update various images to Ver.5.2.6.</li> </ul>  |
| F        | 2013.9  | <ul style="list-style-type: none"> <li>• Apply to RDSR.</li> <li>• Update various images to Ver.5.2.11.</li> </ul>  |
| G        | 2013.9  | <ul style="list-style-type: none"> <li>• Add description about sending RDSR.</li> <li>• Update “Special Information” in chapter 2.</li> <li>• Update “Periodic Maintenance” in chapter 14.</li> </ul>   |
| J        | 2014.6  | <ul style="list-style-type: none"> <li>• Apply to SCORE PRO Advance.</li> <li>• Update the labels on DAR.</li> <li>• Apply to RoHS.</li> </ul>  |
| K        | 2015.1  | <ul style="list-style-type: none"> <li>• Update cable and accessory lists.</li> <li>• Update labels.</li> </ul>   |
| K        | 2015.1  | <ul style="list-style-type: none"> <li>• Update cable and accessory lists.</li> <li>• Update labels.</li> </ul>   |

## Revision History

| REVISION | DATE   | COMMENT   |
|----------|--------|---|
| L        | 2015.4 | <ul style="list-style-type: none"> <li>• Apply to Ver.6.0.</li> <li>• Add [Configure Display Unit] dialog in Problem and Action for “Startup the System” in chapter 13.</li> <li>• Add periodic replacement parts.</li> <li>• Add description of pediatrics.</li> <li>• Add Dose Meter Configuration.</li> <li>• Add description on TraceMAP.</li> <li>• Add description of contour enhancement.</li> <li>• Add description of ROI of SIMAP.</li> <li>• Add description of guide display.</li> </ul>  |
| M        | 2015.5 | <ul style="list-style-type: none"> <li>• Add [Geometry] in [Physicians] tab on configuration window.</li> <li>• Add [Geometry] in [Physicians] tab on configuration window.</li> </ul>  |
| N        | 2015.6 | <ul style="list-style-type: none"> <li>• Update [Physicians] and [DICOM] tab on configuration window in chapter 12.</li> </ul>  |
| P        | 2015.8 | <ul style="list-style-type: none"> <li>• Update description of Fluoro MAP, SIMAP (LIVE) and SIMAP (Sub).</li> </ul>   |
| R        | 2016.4 | <ul style="list-style-type: none"> <li>• Apply to Ver.6.2.</li> <li>• Update description of MAP.</li> <li>• Update functions of IVR NEO buttons.</li> <li>• Update Keyboard Shortcuts.</li> <li>• Update the section of Displaying RDSR.</li> <li>• Update description of Bi-plane Loop.</li> <li>• Update the configuration dialog box in chapter 12.</li> <li>• Add description of Processed Image Transfer.</li> <li>• Update the names of buttons and items.</li> <li>• Add [Temporal DSA] for FLUO and RAD DUP setting.</li> <li>• Update the section of Other in chapter 13.</li> </ul> |
| T        | 2016.5 | <ul style="list-style-type: none"> <li>• Update the section of Displaying RDSR.</li> </ul>  |

| REVISION | DATE    | COMMENT   |
|----------|---------|---|
| V        | 2017.5  | <ul style="list-style-type: none"> <li>• Add description of Operator.</li> <li>• Update 4.8.1 IVR NEO Buttons.</li> <li>• Update 12.9 Options Configuration.</li> <li>• Add description of MPPS Support to section 4.2, 4.3.2, 6.2.2, and 12.8.2.</li> <li>• Add 4.11.</li> <li>• Update 13.1 Error Messages.</li> <li>• Add section 13.3 Error Report.</li> <li>• Update 13.4 Other.</li> <li>• Update 12.9.4 Menus and DUP Configuration.</li> </ul>  |
| Y        | 2017.9  | <ul style="list-style-type: none"> <li>• Apply to DAR-9500f S/W Ver.6.5.</li> <li>• Update main window.</li> <li>• Update DSA mode in section 2.5.1.</li> <li>• Update cable and accessory lists.</li> <li>• Update labels.</li> <li>• Add Warning and Caution for Precautions in Use.</li> <li>• Add description of Wireless Communication in chapter2.</li> <li>• Add description of Auto Stitching in chapter 2 and section 4.4.6.</li> <li>• Add description of SMART Touch in section 14.4.</li> <li>• Update New Study windows in section 4.3.</li> <li>• Update series tab of Studies Management window.</li> <li>• Update a dialog box of Media Write.</li> <li>• Update DSA Tool bar.</li> <li>• Update [Hardware] tab and [Physicians]-[Geometry Configuration] tab of option configuration window.</li> <li>• Add 4.11 SMART Touch (Either-or IVR NEO/IVR Shuttle).</li> </ul> |
| AA       | 2017.11 | <ul style="list-style-type: none"> <li>• Add SMART Touch (Option) list in section 1.5.</li> <li>• Add standard components of SMART Touch in section 2.3.1.</li> <li>• Add description of Option button and Change of Preset in section 4.11.1.</li> <li>• Add [Touch] tab on [Physician] tab of option configuration window.</li> <li>• Add description of Flex-APS (Option) in section 2.5.3, 4.4.6, 4.8.1, 4.11.6 and 9.3.</li> <li>• Add description of Error Messages Related to SMART Touch in section 13.1</li> <li>• Add 13.5 Displaying Image of SMART Display (Option).</li> </ul>   |

## Revision History

| REVISION | DATE    | COMMENT  |
|----------|---------|--|
| AB       | 2018.1  | <ul style="list-style-type: none"> <li>Add Chapter 11 SCORE StentView (Option).</li> <li>Add Chapter 12 SCORE StentShot (Option).</li> <li>Add Chapter 13 Temporary License.</li> </ul>  |
| AC       | 2018.4  | <ul style="list-style-type: none"> <li>Add "Note" in section 4.11.4 Monitor Panel (Only when SMART Display is connected).</li> <li>Revise "Note" for Auto Stitching (Option) in section 4.4.6.</li> </ul>  |
| AD       | 2018.8  | <ul style="list-style-type: none"> <li>Apply to DAR-9500f S/W Rev.6.7.</li> <li>Add Chapter 13 SCORE Chase (Option).</li> <li>Add Chapter 14 SMART Display (Option).</li> <li>Add section 13.4.5 SPOT-DSA and 13.4.8 Change of Display.</li> <li>Add description of SPOT-DSA in section 13.1.2, 13.3.1, 13.3.2.1, and 13.4.10.</li> <li>Add "Note" in section 13.4.3 and 13.4.4.</li> </ul>  |
| AE       | 2018.11 | <ul style="list-style-type: none"> <li>Add "Warning" and "Caution" in Precautions in Use.</li> <li>Add "Note" in section 2.6.4.</li> </ul>   |
| AF       | 2018.11 | <ul style="list-style-type: none"> <li>Modify the image defects.</li> </ul>  |
| AG       | 2018.12 | <ul style="list-style-type: none"> <li>Apply to DAR-9500f S/W Rev.6.8.</li> <li>Add "Note" of X-ray generator activation in section 3.1 and 4.3.2.</li> <li>Add description of "Show dose warning" in section 4.2 and 4.4.1.</li> <li>Change Set Acquisition Parameters window in section 4.4.8.</li> <li>Modify default value of Brightness and Contrast in section 7.2.1 and 7.2.2.</li> <li>Change Quantitative Coronary Analysis window, and add procedures from step 16 to 19 in section 8.2.</li> <li>Update [Notification] tab and add description of dose in section 17.8.3.</li> <li>Update [Devices] tab and add description of Screen Saver in section 17.8.3.</li> </ul> |
| AH       | 2019.1  | <ul style="list-style-type: none"> <li>Add "Note" in section 2.6.8.</li> <li>Revise "Note" in section 3.1.</li> <li>Add periodic replacement parts in section 19.3.2.</li> </ul>   |
| AJ       | 2019.1  | <ul style="list-style-type: none"> <li>Revise "Note" in section 2.6.8.</li> <li>Delete sentences in section 3.1.</li> <li>Modify the list in section 4.2.</li> </ul>   |



| REVISION | DATE       | COMMENT   |
|----------|------------|---|
| AK       | 2019.5     | <ul style="list-style-type: none"> <li>• Update the cable list.</li> <li>• Delete labels in section 2.6.2.</li> <li>• Update [DICOM] and [Hardware] tabs in section 17.9.3.</li> <li>• Update [Hardware] tab in section 10.5.</li> </ul>  |
| AL       | 2019.8     | <ul style="list-style-type: none"> <li>• Apply to DAR-9500f S/W Rev.6.9.</li> </ul>   |
| AM       | 2020.7.31  | <ul style="list-style-type: none"> <li>• Modify periodic replacement parts in section 19.3.2.</li> </ul>  |
| AN       | 2021.2.8   | <ul style="list-style-type: none"> <li>• Add an item to "Image Pre-Processing" in the table of DUP parameters in section 17.9.4.</li> </ul>   |
| AP       | 2022.11.30 | <ul style="list-style-type: none"> <li>• Add a rating label in section 2.6.2.</li> </ul>  |
| AR       | 2023.5.8   | <ul style="list-style-type: none"> <li>• Update section 4.4.9 Manually Resetting the Dose.</li> </ul>   |
| AT       | 2024.4.30  | <ul style="list-style-type: none"> <li>• Add section 3.2 User Authentication.</li> <li>• Add section 4.5.11 Circle Display.</li> <li>• Add procedure in section 9.2.</li> <li>• Update [System] tab in section 10.5.</li> <li>• Add "NOTE" in section 13.4.6.</li> <li>• Update [Display] tab in section 17.8.3.</li> <li>• Add [Cyber Security] tab in section 17.8.3.</li> <li>• Add section 18.4 Unable to Sign In.</li> </ul> |

---

This page is intentionally left blank.

# Table of Contents

|  |      |
|--|------|
| Introduction .....                                   | i    |
| Precautions in Use .....                             | iv   |
| Connector .....                                      | vii  |
| Other Software Installation .....                    | vii  |
| Internet Connection .....                            | vii  |
| Scan Converter BNC Cable Positioning .....           | viii |
| Data Erasure or Loss .....                           | viii |
| Disposal .....                                       | viii |
| Cautions on Wireless Communication Using Mouse ..... | ix   |
| Network Connection .....                             | x    |
| Action for Environment .....                         | xi   |
| Pediatric Imaging .....                              | xii  |
| General .....  | xii  |
| Suggestions for Pediatric Imaging .....              | xii  |
| The Limited Product Warranty .....                   | xiv  |
| Revision History .....                               | xv   |

## Chapter 1 About this Manual

---

|                               |     |
|-------------------------------|-----|
| 1.1 Safety .....              | 1-2 |
| 1.2 Manual Organization ..... | 1-2 |
| 1.3 Operation Modes .....     | 1-2 |
| 1.4 Operator Profile .....    | 1-3 |
| 1.5 Terminology .....         | 1-4 |
| FPD .....                     | 1-4 |
| SMART Display (Option) .....  | 1-5 |

## Chapter 2 System Overview

---

|                                 |     |
|---------------------------------|-----|
| 2.1 Introduction .....          | 2-2 |
| 2.2 Features .....              | 2-6 |
| 2.3 System Components .....     | 2-7 |
| 2.3.1 Standard Components ..... | 2-7 |
| IVR NEO .....                   | 2-7 |

|       |  |             |
|-------|--|-------------|
|       | SMART Touch  | 2-7         |
| 2.3.2 | Optional Components  | 2-8         |
| 2.4   | Specifications   | 2-9         |
| 2.4.1 | Image Processing   | 2-9         |
| 2.4.2 | Image Display  | 2-9         |
| 2.4.3 | Image Recording  | 2-10        |
| 2.4.4 | Network  | 2-10        |
| 2.5   | Function   | 2-11        |
| 2.5.1 | Image Acquisition  | 2-11        |
| 2.5.2 | Fluoroscopy  | 2-12        |
| 2.5.3 | Image Processing   | 2-13        |
| 2.5.4 | Image Analysis   | 2-13        |
| 2.5.5 | Network  | 2-14        |
| 2.5.6 | Image Recording  | 2-14        |
| 2.5.7 | C-arm Reposition   | 2-14        |
| 2.5.8 | Table Reposition   | 2-14        |
| 2.6   | Special Information  | 2-15        |
| 2.6.1 | Classification   | 2-15        |
|       | The classification of protection against electric shock:   | 2-15        |
|       | The degree of protection against electric shock:   | 2-15        |
|       | The classification of <b>EMC (Electro-Magnetic Compatibility)</b> :  | <b>2-15</b> |
|       | The degree of protection against effects of water:   | 2-16        |
|       | The degree of safety of application in the presence of a flammable<br>anesthetic mixture with air or with oxygen or nitrous oxide: | 2-16        |
| 2.6.2 | Labeling   | 2-17        |
|       | Warning/Caution Labels   | 2-17        |
|       | Storage and Transport Labels   | 2-19        |
|       | Others   | 2-19        |
| 2.6.3 | Symbols  | 2-23        |
| 2.6.4 | Environmental Conditions   | 2-24        |
|       | Setup and Operating  | 2-24        |
|       | Storage and Transportation   | 2-24        |
|       | Conformity of Equipment  | 2-25        |
| 2.6.5 | EMC (Electro Magnetic Compatibility) Information   | 2-25        |
|       | Operating Condition  | 2-25        |
| 2.6.6 | Wireless Communications  | 2-25        |
|       | Wireless Mouse   | 2-25        |
|       | Cable List   | 2-26        |
|       | Guidance and Manufacturer's Declaration-Electromagnetic Emissions  | 2-39        |
|       | Guidance and Manufacturer's Declaration-Electromagnetic Immunity   | 2-40        |
|       | Guidance and Manufacturer's Declaration-Electromagnetic Immunity   |             |

|  |      |
|--|------|
|  | 2-42 |
| Recommended Separation Distance Between Portable and Mobile RF Communication Equipment and the EQUIPMENT or SYSTEM . | 2-44 |
| Accessory List . . . . .   | 2-46 |
| Essential Performance . . . . .  | 2-47 |
| 2.6.7 Statement of Compliance (for Europe) . . . . .   | 2-48 |
| Regulatory Information . . . . .   | 2-48 |
| Company's Quality System . . . . .   | 2-48 |
| 2.6.8 Other . . . . .  | 2-48 |

## Chapter 3 System Startup and Shutdown

---

|  |      |
|--|------|
| 3.1 Startup . . . . .                      | 3-2  |
| 3.2 User Authentication . . . . .          | 3-4  |
| 3.2.1 User Level . . . . .                 | 3-4  |
| 3.2.2 Sign In . . . . .                    | 3-4  |
| 3.2.3 Input User ID and Password . . . . . | 3-5  |
| 3.2.4 Change Password . . . . .            | 3-7  |
| 3.2.5 Sign Out . . . . .                   | 3-8  |
| 3.2.6 Emergency Sign In . . . . .          | 3-8  |
| 3.3 Shutdown . . . . .                     | 3-10 |
| 3.4 Power Failure . . . . .                | 3-13 |

## Chapter 4 Image Acquisition

---

|  |      |
|--|------|
| 4.1 Introduction . . . . .                                 | 4-2  |
| 4.2 The Main User Interface Window . . . . .               | 4-2  |
| 4.2.1 Reference and Acquisition Monitors . . . . .         | 4-8  |
| 4.3 Entering Study Information . . . . .                   | 4-10 |
| 4.3.1 Introduction . . . . .                               | 4-10 |
| 4.3.2 Defining a New Study . . . . .                       | 4-10 |
| Define New Study with Worklist Server Management . . . . . | 4-11 |
| Define New Study with Local Worklist Management . . . . .  | 4-15 |
| Define New Study with Direct Entry . . . . .               | 4-17 |
| Patient Orientation . . . . .                              | 4-19 |
| Display New Study Information . . . . .                    | 4-21 |
| 4.4 Making Acquisition . . . . .                           | 4-23 |

|        |   |      |
|--------|---|------|
| 4.4.1  | Introduction  | 4-23 |
| 4.4.2  | Adding a Comment  | 4-25 |
| 4.4.3  | Choosing Acquisition Modes                              | 4-26 |
| 4.4.4  | Making Live Fluoro Acquisitions                         | 4-30 |
|        | FluoroMAP   | 4-31 |
|        | DSA-MAP   | 4-32 |
|        | BlankMAP  | 4-34 |
|        | LIVE Mode   | 4-35 |
|        | TraceMAP  | 4-37 |
| 4.4.5  | Recording Fluoro  | 4-38 |
|        | Direct Fluoro Record                                    | 4-38 |
|        | Last N Seconds Fluoro Record                            | 4-39 |
| 4.4.6  | Recording Rad   | 4-40 |
|        | DA  | 4-40 |
|        | DSA   | 4-41 |
|        | Flex-APS (Option)                                       | 4-41 |
|        | RSM-DSA   | 4-42 |
|        | ROT-DA and ROT-RSM-DSA                                  | 4-42 |
|        | ROT-DSA   | 4-43 |
|        | SPOT  | 4-45 |
|        | Precession/pendulum DA and RSM-DSA                      | 4-45 |
|        | Stage Acquisition                                       | 4-46 |
| 4.4.7  | Flipping Images   | 4-47 |
| 4.4.8  | Changing the Acquisition Parameter Settings             | 4-47 |
| 4.4.9  | Manually Resetting the Dose                             | 4-49 |
| 4.5    | Actions When Study Open for Acquisition                 | 4-50 |
| 4.5.1  | Continuing a Study                                      | 4-50 |
| 4.5.2  | Using Select  | 4-51 |
| 4.5.3  | Using Split   | 4-52 |
| 4.5.4  | Filtering Images by C-arm Position                      | 4-53 |
| 4.5.5  | Moving C-arm to Match Selected Image (Reposition)       | 4-54 |
| 4.5.6  | Moving Table to Match Selected Image (Table Reposition) | 4-55 |
| 4.5.7  | Sending Selected Image                                  | 4-55 |
| 4.5.8  | Deleting Selected Image                                 | 4-56 |
| 4.5.9  | Displaying Previous Studies                             | 4-56 |
| 4.5.10 | Sketch Display  | 4-57 |
| 4.5.11 | Circle Display  | 4-59 |
| 4.6    | Working With Reference Images                           | 4-60 |
| 4.6.1  | Introduction  | 4-60 |
| 4.6.2  | Adding a Reference Image                                | 4-60 |
| 4.6.3  | Displaying Images Selected as Reference                 | 4-61 |
| 4.6.4  | Deleting a Reference Image                              | 4-61 |
| 4.6.5  | Filtering a Reference Image                             | 4-62 |

|        |   |       |
|--------|---|-------|
| 4.7    | Acquisition Keys and GUI Buttons                              | 4-63  |
|        | Menu Interaction (If Digital Program Menu is not displayed)   | 4-63  |
|        | Menu Interaction (If Rad, Fluoro, Comment Menu are displayed) | 4-64  |
|        | Reference Images  | 4-64  |
|        | 16 Up Select  | 4-64  |
| 4.8    | IVR NEO (Either-or SMART Touch)                               | 4-65  |
| 4.8.1  | IVR NEO Buttons   | 4-66  |
| 4.8.2  | IVR NEO Joystick  | 4-71  |
| 4.9    | Keyboard ShortCuts  | 4-72  |
| 4.10   | IVR Shuttle (Option) (Either-or SMART Touch)                  | 4-73  |
| 4.10.1 | Functions of IVR Shuttle Buttons and Controls                 | 4-74  |
| 4.11   | SMART Touch (Either-or IVR NEO/IVR Shuttle)                   | 4-77  |
| 4.11.1 | Custom Panel  | 4-79  |
|        | Option Button   | 4-79  |
|        | Modify Assignment of Buttons                                  | 4-84  |
|        | Change of Button Colors                                       | 4-88  |
| 4.11.2 | Program Panel   | 4-90  |
|        | Switch RAD DUP  | 4-90  |
|        | Switch RAD Plane (In case of Bi-plane system)                 | 4-90  |
|        | Switch Fuoro DUP  | 4-90  |
|        | Switch FOV  | 4-90  |
| 4.11.3 | Tool Panel  | 4-91  |
| 4.11.4 | Monitor Panel (Only when SMART Display is connected)          | 4-93  |
|        | Change of Preset  | 4-94  |
|        | Change of Segment   | 4-94  |
|        | Change of Video Input   | 4-96  |
| 4.11.5 | Side Menu   | 4-98  |
|        | X-ray OFF   | 4-98  |
|        | Function Button   | 4-99  |
|        | Stop Fluoro Alarm   | 4-99  |
| 4.11.6 | SMART Touch Buttons   | 4-100 |
| 4.11.7 | SMART Touch Joystick  | 4-104 |
| 4.12   | Closing the Active Study                                      | 4-105 |
| 4.13   | Dose Report Image   | 4-107 |

## Chapter 5 Review Basics

---

|     |                             |     |
|-----|-----------------------------|-----|
| 5.1 | Startup and Study Selection | 5-2 |
| 5.2 | Viewing Images              | 5-3 |

|       |                                    |     |
|-------|------------------------------------|-----|
| 5.2.1 | Controlling Image Playback         | 5-5 |
| 5.2.2 | Adjusting Image Appearance         | 5-5 |
| 5.2.3 | Identifying Loops and Still Images | 5-5 |
| 5.3   | Displaying Study Information       | 5-6 |
| 5.4   | Creating Image Annotation          | 5-6 |
| 5.5   | Saving an Image as a New Object    | 5-7 |
| 5.6   | Printing Images                    | 5-7 |
| 5.7   | Using Advanced Features            | 5-8 |
| 5.8   | Completing Your Review             | 5-8 |

## Chapter 6 Studies Management

---

|       |   |      |
|-------|---|------|
| 6.1   | The Studies Management Window           | 6-2  |
| 6.1.1 | Study Location                          | 6-3  |
| 6.1.2 | Action Bar                              | 6-4  |
| 6.1.3 | Search Criteria                         | 6-5  |
| 6.1.4 | Study List                              | 6-5  |
|       | Sorting                                 | 6-6  |
|       | Scrolling in the Study List             | 6-6  |
|       | Widen or Narrow Columns as Desired      | 6-6  |
|       | Making Selections in the Study List     | 6-6  |
| 6.1.5 | Series Tab                              | 6-7  |
| 6.1.6 | Options Tab                             | 6-7  |
| 6.1.7 | Messages & Information                  | 6-8  |
| 6.2   | Managing Studies                        | 6-9  |
| 6.2.1 | Introduction                            | 6-9  |
| 6.2.2 | Finding Studies                         | 6-10 |
| 6.2.3 | Downloading Studies from Servers        | 6-12 |
| 6.2.4 | Displaying RDSR                         | 6-13 |
| 6.2.5 | Importing Studies from CD/DVD           | 6-14 |
| 6.2.6 | Sending Studies                         | 6-14 |
| 6.2.7 | Writing to CD/DVD                       | 6-15 |
|       | Recommended Media (As of January, 2015) | 6-15 |
|       | Select the Studies                      | 6-16 |
|       | Optionally Select Fewer Items           | 6-17 |
|       | Perform the Write                       | 6-18 |
|       | Insert Additional Blank CD/DVDs         | 6-19 |
|       | Cancel a Media Write Operation          | 6-19 |
|       | Label CD/DVD                            | 6-19 |



|       |                                |      |
|-------|--------------------------------|------|
| 6.2.8 | Setting Study Protection ..... | 6-19 |
| 6.2.9 | Deleting Local Studies .....   | 6-20 |

## Chapter 7 Review Reference

---

|       |  |      |
|-------|--|------|
| 7.1   | Controlling Image Playback .....             | 7-2  |
| 7.1.1 | Mouse and Keyboard .....                     | 7-2  |
|       | Play Acquisition Monitor with Keyboard ..... | 7-2  |
| 7.1.2 | Loop Playback Mode .....                     | 7-3  |
| 7.1.3 | Adjusting Playback Speed .....               | 7-3  |
| 7.2   | Adjusting Image Appearance .....             | 7-4  |
| 7.2.1 | Basic Brightness/Contrast Control .....      | 7-4  |
|       | Using Side Menu Buttons .....                | 7-4  |
|       | Using the Mouse .....                        | 7-4  |
| 7.2.2 | Advanced Brightness/Contrast Control .....   | 7-5  |
|       | Auto Window Level Feature .....              | 7-5  |
|       | LUT Feature .....                            | 7-6  |
| 7.2.3 | Applying Filters .....                       | 7-7  |
| 7.2.4 | Zooming and Panning .....                    | 7-8  |
|       | Using the Zoom Button and Mouse .....        | 7-8  |
| 7.3   | Working with Multiple Studies .....          | 7-9  |
| 7.4   | Displaying Study Information .....           | 7-10 |
| 7.5   | Modifying Study Information .....            | 7-12 |
| 7.6   | Anonymizing Study Information .....          | 7-15 |
| 7.7   | Working with Image Annotations .....         | 7-17 |
| 7.7.1 | Measurement Tools .....                      | 7-19 |
|       | Caliper .....                                | 7-19 |
|       | Line and Curve .....                         | 7-19 |
| 7.8   | Saving Images .....                          | 7-20 |
| 7.9   | Print Images .....                           | 7-21 |
| 7.9.1 | Selecting Images for Printing .....          | 7-21 |
| 7.9.2 | Printing Basics .....                        | 7-22 |
| 7.9.3 | Print Features .....                         | 7-23 |
|       | Horizontal/Vertical .....                    | 7-23 |
|       | Modifying the Header and Footer .....        | 7-23 |
|       | Modifying the Image Legend .....             | 7-29 |
|       | Selecting Page Range .....                   | 7-29 |
|       | Setting Printer LUT .....                    | 7-30 |
|       | Setting Printer .....                        | 7-30 |

|       |                                |      |
|-------|--------------------------------|------|
| 7.9.4 | Setting DICOM Printer .....    | 7-31 |
| 7.10  | Modifying Image Comments ..... | 7-33 |

## Chapter 8 Performing Angiographic Analysis

---

|       |                                      |      |
|-------|--------------------------------------|------|
| 8.1   | Introduction .....                   | 8-2  |
| 8.2   | Quantitative Coronary Analysis ..... | 8-2  |
| 8.3   | Left Ventricular Analysis .....      | 8-7  |
| 8.4   | Image Calibration Methods .....      | 8-11 |
| 8.4.1 | Using ISO Center Value .....         | 8-11 |
| 8.4.2 | Using a Catheter (QCA) .....         | 8-11 |
| 8.4.3 | Using Reference Distance (QCA) ..... | 8-12 |
| 8.4.4 | Using a Grid (QCA and LV) .....      | 8-13 |
| 8.4.5 | Using a Reference Ball (LV) .....    | 8-13 |

## Chapter 9 DSA Image Adjustment Tools

---

|     |  |      |
|-----|--|------|
| 9.1 | Introduction .....                     | 9-2  |
| 9.2 | Adjusting Subtraction Parameters ..... | 9-3  |
| 9.3 | Reregistration .....                   | 9-5  |
| 9.4 | Peak Hold .....                        | 9-6  |
| 9.5 | Contour Enhancement .....              | 9-8  |
| 9.6 | Mask Region .....                      | 9-11 |
| 9.7 | Adjusting Shutter .....                | 9-13 |

## Chapter 10 Difference with Bi-plane Option

---

|        |  |      |
|--------|--|------|
| 10.1   | Introduction .....                         | 10-2 |
| 10.1.1 | Standard Components .....                  | 10-2 |
| 10.1.2 | Bi-plane User Interface .....              | 10-2 |
| 10.1.3 | Reference Monitors .....                   | 10-4 |
| 10.2   | Bi-plane Image Acquisition .....           | 10-5 |
| 10.2.1 | Entering Study Information .....           | 10-5 |
| 10.2.2 | Acquiring Loops when Bi-plane Active ..... | 10-5 |

|  |       |
|--|-------|
| 10.2.3 Bi-plane Loop .....                               | 10-7  |
| Reference Monitor .....                                  | 10-7  |
| Acquisition Monitor .....                                | 10-8  |
| 10.2.4 Filtering Bi-plane Images by C-arm Position ..... | 10-12 |
| 10.3 Bi-plane Image Preview .....                        | 10-14 |
| 10.4 Shutdown .....                                      | 10-15 |
| 10.5 Administration .....                                | 10-16 |

## Chapter 11 SCORE StentView (Option)

---

|   |       |
|---|-------|
| 11.1 Overview .....                           | 11-2  |
| 11.1.1 Features .....                         | 11-2  |
| 11.2 Specification .....                      | 11-3  |
| 11.2.1 Specification of SCORE StentView ..... | 11-3  |
| 11.2.2 Features .....                         | 11-3  |
| Image Acquisition .....                       | 11-3  |
| Image Processing .....                        | 11-3  |
| 11.3 Handling Method .....                    | 11-4  |
| 11.3.1 Introduction .....                     | 11-4  |
| 11.3.2 Displaying StentView .....             | 11-4  |
| Full Screen View .....                        | 11-4  |
| 11.3.3 StentView Processing .....             | 11-4  |
| Setting ROI .....                             | 11-7  |
| 11.3.4 Change View .....                      | 11-11 |
| 11.3.5 Setting StentView .....                | 11-11 |
| Setting Shortcut Keys .....                   | 11-12 |
| Setting StentView .....                       | 11-14 |
| Geometry .....                                | 11-16 |
| 11.4 Troubleshooting .....                    | 11-18 |
| 11.4.1 Phenomenon and Action .....            | 11-18 |

## Chapter 12 SCORE StentShot (Option)

---

|   |      |
|---|------|
| 12.1 Overview .....                           | 12-2 |
| 12.1.1 Features .....                         | 12-2 |
| 12.2 Specification .....                      | 12-3 |
| 12.2.1 Specification of SCORE StentShot ..... | 12-3 |

|                                    |       |
|------------------------------------|-------|
| 12.2.2 Features .....              | 12-3  |
| Image Acquisition .....            | 12-3  |
| Image Processing .....             | 12-3  |
| 12.3 Handling Method .....         | 12-4  |
| 12.3.1 Introduction .....          | 12-4  |
| 12.3.2 Displaying StentShot .....  | 12-4  |
| Full Screen View .....             | 12-4  |
| 12.3.3 StentShot Processing .....  | 12-4  |
| Setting ROI .....                  | 12-7  |
| 12.3.4 Change View .....           | 12-11 |
| 12.3.5 Setting StentShot .....     | 12-11 |
| Setting Shortcut Keys .....        | 12-12 |
| Setting StentShot .....            | 12-14 |
| Geometry .....                     | 12-16 |
| 12.4 Troubleshooting .....         | 12-18 |
| 12.4.1 Phenomenon and Action ..... | 12-18 |

## Chapter 13 SCORE Chase (Option)

---

|  |       |
|--|-------|
| 13.1 Overview .....  | 13-2  |
| 13.1.1 Purpose .....   | 13-2  |
| 13.1.2 Feature .....   | 13-2  |
| 13.2 Components .....  | 13-3  |
| 13.2.1 Chase Console .....                                   | 13-3  |
| 13.3 Specification .....                                     | 13-5  |
| 13.3.1 Specification of SCORE Chase .....                    | 13-5  |
| 13.3.2 Features .....  | 13-6  |
| Image Acquisition .....                                      | 13-6  |
| Image Processing .....                                       | 13-6  |
| 13.4 Handling Method .....                                   | 13-7  |
| 13.4.1 Introduction .....                                    | 13-7  |
| 13.4.2 Manual Chase .....                                    | 13-7  |
| 13.4.3 RemoteChase .....                                     | 13-9  |
| 13.4.4 RemoteChase-DSA .....                                 | 13-12 |
| 13.4.5 SPOT-DSA .....  | 13-15 |
| 13.4.6 PostProcess Stitch .....                              | 13-17 |
| 13.4.7 Positioning .....                                     | 13-19 |
| 13.4.8 Change of Display .....                               | 13-21 |
| 13.4.9 Image Loop and Identification of Stitched Image ..... | 13-22 |
| 13.4.10 Setting of SCORE Chase .....                         | 13-23 |
| Default Setting of Table Position Display .....              | 13-23 |

|   |       |
|---|-------|
| Setting IVR NEO and Shortcut Keys ..... | 13-24 |
| Setting of SMART Touch .....            | 13-25 |

## Chapter 14 SMART Display

---

|                                     |       |
|-------------------------------------|-------|
| 14.1 Introduction .....             | 14-2  |
| 14.2 Handling Method .....          | 14-4  |
| 14.2.1 Main Screen .....            | 14-4  |
| 14.2.2 Basic Settings .....         | 14-5  |
| 14.2.3 Preset .....                 | 14-6  |
| 14.2.4 Video Input .....            | 14-7  |
| 14.2.5 Switch Segment .....         | 14-7  |
| 14.2.6 Switch Video Input .....     | 14-8  |
| 14.2.7 Creating Preset .....        | 14-9  |
| 14.3 Setting of SMART Display ..... | 14-12 |
| Setting of IP Address .....         | 14-12 |
| 14.4 Trouble Shooting .....         | 14-13 |
| 14.4.1 Error Message .....          | 14-13 |
| 14.4.2 Displaying Image .....       | 14-14 |

## Chapter 15 Temporary License

---

|                                      |      |
|--------------------------------------|------|
| 15.1 Introduction .....              | 15-2 |
| 15.2 Start a Temporary License ..... | 15-3 |
| 15.3 End a Temporary License .....   | 15-4 |

## Chapter 16 Mini Viewer

---

|                            |      |
|----------------------------|------|
| 16.1 The Mini Viewer ..... | 16-2 |
| Install DirectX 11 .....   | 16-4 |

## Chapter 17 Administration

---

|        |   |       |
|--------|---|-------|
| 17.1   | Introduction                              | 17-2  |
| 17.2   | Switching Modes                           | 17-2  |
| 17.2.1 | Switch to Administration Mode             | 17-2  |
| 17.2.2 | Switch to User Mode                       | 17-2  |
| 17.3   | Viewing User Type                         | 17-3  |
| 17.4   | Storage Commit Messages                   | 17-4  |
| 17.5   | Working with LUT Profiles                 | 17-4  |
| 17.5.1 | Introduction                              | 17-4  |
| 17.5.2 | LUT Profile Definition                    | 17-5  |
| 17.5.3 | Assigning Display Adjustment LUT Profiles | 17-7  |
| 17.5.4 | Assigning Printer LUT Profiles            | 17-8  |
| 17.5.5 | Setting the Default Printer LUT Profile   | 17-9  |
| 17.6   | DICOM Monitor Linearization               | 17-10 |
| 17.7   | Options Configuration                     | 17-14 |
| 17.7.1 | Introduction                              | 17-14 |
| 17.7.2 | Options Tab                               | 17-14 |
| 17.7.3 | Configuration Dialog Box                  | 17-15 |
|        | System Tab                                | 17-17 |
|        | DICOM Hosts Tab                           | 17-20 |
|        | DICOM Tab                                 | 17-24 |
|        | Software Tab                              | 17-27 |
|        | Hardware Tab                              | 17-30 |
|        | Hardware (Shutter) Tab                    | 17-32 |
|        | Storage (Local) Tab                       | 17-33 |
|        | Storage (Network) Tab                     | 17-35 |
|        | Display Tab                               | 17-37 |
|        | Menus Tab                                 | 17-41 |
|        | Notifications Tab                         | 17-42 |
|        | Physicians Tab                            | 17-44 |
|        | Study Information Tab                     | 17-50 |
|        | Fonts and Colors Tab                      | 17-52 |
|        | Database Tab                              | 17-54 |
|        | Devices Tab                               | 17-56 |
|        | External Software Tab                     | 17-58 |
|        | Fusion Tab                                | 17-59 |
|        | Logs Tab                                  | 17-59 |
|        | Cyber Security Tab                        | 17-60 |
| 17.7.4 | Menus and DUP Configuration               | 17-69 |
|        | Configure the Comment Menu                | 17-70 |
|        | Configure Fluoro/Rad Menus and DUPs       | 17-71 |
|        | Configure Fluoro/Rad Menus                | 17-71 |

|                                 |       |
|---------------------------------|-------|
| Configure Fluoro/Rad DUPs ..... | 17-73 |
|---------------------------------|-------|

## Chapter 18 Troubleshooting

---

|   |       |
|---|-------|
| 18.1 Error Messages .....                   | 18-2  |
| Error Messages Related to FPD .....         | 18-3  |
| Error Messages Related to SMART Touch ..... | 18-3  |
| 18.2 Startup the System .....               | 18-5  |
| 18.3 Error Report .....                     | 18-7  |
| 18.4 Unable to Sign In .....                | 18-10 |
| 18.5 Other .....                            | 18-11 |

## Chapter 19 Maintenance

---

|                                     |      |
|-------------------------------------|------|
| 19.1 Introduction .....             | 19-2 |
| 19.2 Daily Maintenance .....        | 19-3 |
| 19.2.1 Startup Maintenance .....    | 19-3 |
| Purpose .....                       | 19-3 |
| Procedures .....                    | 19-3 |
| 19.2.2 Finishing Operations .....   | 19-4 |
| Post-operation Maintenance .....    | 19-4 |
| 19.3 Periodic Maintenance .....     | 19-5 |
| 19.3.1 Calibration .....            | 19-5 |
| FPD Calibration .....               | 19-5 |
| 19.3.2 Preventive Maintenance ..... | 19-6 |
| Periodic Part Replacement .....     | 19-6 |

## Chapter 20 Appendix

---

|                                 |      |
|---------------------------------|------|
| 20.1 DICOM Concepts .....       | 20-2 |
| 20.1.1 The DICOM Standard ..... | 20-2 |
| 20.2 Glossary .....             | 20-3 |

This page is intentionally left blank.



# *Chapter 1*

---

## About this Manual

This chapter introduces this Manual, describing its organization and terminology.

### Description of Chapter



|     |                               |     |
|-----|-------------------------------|-----|
| 1.1 | Safety . . . . .              | 1-2 |
| 1.2 | Manual Organization . . . . . | 1-2 |
| 1.3 | Operation Modes . . . . .     | 1-2 |
| 1.4 | Operator Profile . . . . .    | 1-3 |
| 1.5 | Terminology . . . . .         | 1-4 |

## 1.1 Safety

Read all of the instructions contained in this manual before operating this product and save this manual for later use. Follow all warning, cautions, and instructions.

## 1.2 Manual Organization

This manual is organized in logical sequence that reflects typical product use. Product walkthroughs and user interface descriptions are provided first, followed by reference material and advanced feature descriptions in later chapters. Administration-related material is presented last.

See  ["17 Administration"](#). This manual is intended to be read sequentially through the Studies Management chapter. See  ["6 Studies Management"](#). All later chapters can be referenced on an as-needed basis.

## 1.3 Operation Modes

DAR-9500f features two operation modes; User and Administration. Except when otherwise stated, the information in this manual is for the User.

| Mode           | Description   |
|----------------|---|
| User           | Any user not defined as Administrator runs DAR-9500f in User mode. This is the normal operation in which all day-to-day work should be done. There are usually at least two user types; Technician and Physician. The Technician prepares for a procedure by entering patient and procedure information whereas the Physician performs the procedure. |
| Administration | The system administrator runs DAR-9500f in Administration mode, gaining the ability to search specific servers and CD/DVD drives and perform full configuration.  |

# 1.4 Operator Profile

1

| Item                   | Details   |
|------------------------|---|
| Age                    | Age that person can obtain the license of Radiological Technologist or a license equal to it. |
| Sex                    | No limitation.  |
| Nationality            | No limitation.  |
| Education              | Radiological Technologist or person who has a license equal to it.                            |
| Knowledge              | Radiological Technologist or person who has a license equal to it.                            |
| Language               | Can read and understand English.  |
| Experience             | Every operator needs to take training for operating the equipment before using the equipment. |
| Permissible Impairment | Corrected visibility is over 0.7 in the decimal number.                                       |

# 1.5 Terminology

The following terminology is used in this manual:

| Terms                                 | Description   |
|---------------------------------------|---|
| GUI                                   | Graphical User Interface, a graphical system of windows, buttons, and controls used to simply communicate with and operate systems.   |
| Image or Loop                         | Used interchangeably, any still image or image sequence acquired during a study. A sequence of images can be played back at acquired speed (e.g., 30fps), providing high-resolution full-motion video.  |
| Study                                 | Refers to both the medical diagnostic procedure itself and all related information and images captured during the procedure. The study also includes descriptive information such as patient demographics and identification codes entered before the procedure is performed. It may also include annotated images. |
| DICOM                                 | Digital Imaging and Communication in Medicine, is an set of rules that allow medical images and associated information to be exchanged between imaging equipment, computers, and hospitals.   |
| Digital Program or DUP                | Used interchangeably, digital user programs define the type of acquisition, either Rad or Fluoro, and various acquisition characteristics.  |
| Lossless Compression                  | Digital compression in which data loss never occurs. Original data (image) condition is achieved after decompression.   |
| Lossless Compression                  | Digital compression which is not fully reversible but typically allows images to retain sufficient detail for analysis.   |
| Worklist                              | Also called Modality Worklist, this is a list of studies (entered earlier by administrative personnel) that should be performed by the physician(s).  |
| MWM Server (Modality Worklist Server) | An information system which provides study order information to the modality (DAR-9500f system).  |
| Performing Physician                  | A physician who performs or is directly-involved in the study.  |
| Referring Physician                   | The physician who originally requested the study that must later review the study result.   |
| Procedure Step                        | The DICOM term of study.  |

## ■ FPD

| Terms               | Description    |
|---------------------|----------------|
| SFD-1612AF          | 16x12 inch FPD |
| SFD-1212/SFD-1212AF | 12 inch FPD    |
| SFD-0808/SFD-0808AF | 8 inch FPD     |

**■ SMART Display (Option)**

| Terms       | Description   |
|-------------|---|
| Video Input | Input image from PC, etc.   |
| Segment     | Display video input on the segment.                                   |
| Layout      | Split display into several segments by layout.                        |
| Preset      | Video input is assigned to the segment within layout based on preset. |

 ["20.2 Glossary" P.20-3](#)

This page is intentionally left blank.

# Chapter 2

---

## System Overview

This chapter provides an overview of the system capabilities and features and describes optional functionality.

### ■ Description of Chapter

|     |                           |      |
|-----|---------------------------|------|
| 2.1 | Introduction .....        | 2-2  |
| 2.2 | Features .....            | 2-6  |
| 2.3 | System Components .....   | 2-7  |
| 2.4 | Specifications .....      | 2-9  |
| 2.5 | Function .....            | 2-11 |
| 2.6 | Special Information ..... | 2-15 |

## 2.1 Introduction

DAR-9500f is a digital X-ray angiography system that automatically acquires digital images with FPD (Flat Panel Detector), applies digital processing such as auto window level, and makes the images available for review. Digital Subtraction Angiography (DSA) provides enhanced imaging of static physiology.

DAR-9500f provides a complete digital imaging system for catheterization (cath) labs that produces high-quality digital image loops that may be stored locally or sent to any DICOM 3.0+ archive server across a network.


Providing full connectivity to DICOM3.0+ networks, DAR-9500f supports lossless-compressed image format. DAR-9500f components including cabinet and monitors, are typically installed in the three rooms of a cath lab as follows.

| No. | Room             | Description  |
|-----|------------------|--|
| 1   | Examination Room | Acquisition and Reference monitor, displaying exactly what is shown on the equivalent monitors in the control room. The physician controls X-ray exposure via foot switches or hand-operated switches. The IVR NEO(*) bedside console provides configurable buttons, LED indicators, and a joystick, or the SMART Touch(*) provides configurable buttons and a joystick.<br><br>A FPD (Flat Panel Detector) is installed on their C-arm system.<br>(*) Either Smart Touch or IVR NEO.  |
| 2   | Control Room     | Digital reference cabinet with Acquisition (live) and Reference monitors, Reference computer with keyboard, mouse, CD/DVD writer and SMART Touch or IVR Shuttle (in option)(*). The Acquisition monitor displays all live image loops from the X-ray system. The Reference monitor is used to playback image loops and do other work such as initiating and completing studies, annotating images, and performing analysis.<br><br>The SMART Touch provides configurable buttons and a joystick.<br>(*) Either Smart Touch or IVR Shuttle (in option). |
| 3   | Equipment Room   | The Digital acquisition cabinet with Acquisition computer, C-arm control cabinet, FPD cooling system, X-ray high voltage control cabinet, X-ray tube-cooling unit, and Maintenance computer.   |




- Single-plane System

### <Examination Room>




\*The C-arm system is "C-arm Support MH-300". It is also possible to place with "Ceiling Suspended C-arm Support MH-200S".

### <Control Room>




## 2 System Overview

- Bi-plane System  
    <Examination Room>




\*The C-arm system is "C-arm Support MH-300" and "Ceiling Suspended C-arm Support MH-400"

### <Control Room>



\*For the optional feature of Live Image During Map (LIDM), a LIVE monitor can be added to the Examination room and the Control room.

## &lt;Equipment



\*For Bi-plane system, there are two sets for each equipment except maintenance computer and transformer. See chapter 10 for details.

# 2.2 Features

DAR-9500f has the following features:

- 1 Image loops are acquired at up to 30 frames per second with a resolution of 1024x1024 pixels.
- 2 Clearer fluoroscopy images and sharper radiographic images are provided through noise-reduction processing and edge enhancement.
- 3 The acquired image is stored on the high-speed image disk in real time, and after acquisition is completed, instantaneous review is available. The previous data can be easily searched by patient or other criteria.
- 4 Rad and Fluoro digital program menus make it easy for the operator to quickly set exposure conditions.
- 5 The IVR supporting environment is offered with the IVR NEO or SMART Touch bedside console operation panel.
- 6 Images and loops can be permanently stored in DICOM format on CD/DVD discs, providing both economy and portability.

## 2.3 System Components

DAR-9500f consists of the following standard and optional components. Refer to "Standard Components" for chapter 10.

 ["10.1.1 Standard Components" P.10-2](#)

2

### 2.3.1 Standard Components

#### IVR NEO

| No. | Name  | Installation Location | Quantity |
|-----|---|-----------------------|----------|
| 1   | Digital reference cabinet: with built-in reference computer     | Control room          | 1 set    |
| 2   | Digital acquisition cabinet: with built-in acquisition computer | Equipment room        | 1 set    |
| 3   | FPD cooling unit  | Equipment room        | 1 set    |
| 4   | 19" LCD monitor   | Control room          | 2 sets   |
|     |   | Examination room      | 2 sets   |
| 5   | Keyboard and mouse  | Control room          | 1 set    |
| 6   | Bedside console IVR NEO   | Examination room      | 1 set    |
| 7   | FPD (Flat Panel Detector)                                       | Examination room      | 1 set    |

#### SMART Touch

| No. | Name  | Installation Location | Quantity |
|-----|---|-----------------------|----------|
| 1   | Digital reference cabinet: with built-in reference computer     | Control room          | 1 set    |
| 2   | Digital acquisition cabinet: with built-in acquisition computer | Equipment room        | 1 set    |
| 3   | FPD cooling unit  | Equipment room        | 1 set    |
| 4   | 19" LCD monitor   | Control room          | 2 sets   |
|     |   | Examination room      | 2 sets   |
| 5   | Keyboard and mouse  | Control room          | 1 set    |
| 6   | SMART Touch   | Examination room      | 1 set    |
| 7   | FPD (Flat Panel Detector)                                       | Examination room      | 1 set    |



The system includes 40GB of memory and 240GB or more of hard drive space.

### 2.3.2 Optional Components

---

| No. | Name                                     | Quantity |
|-----|--|----------|
| 1   | Scan converter                           | 1 set    |
| 2   | Dosimeter                                | 1 set    |
| 3   | Live Image During Map (LIDM)             | 1 set    |
| 4   | IVR Shuttle (Combined with IVR NEO only) | 1 set    |
| 5   | SMART Display                            | 1 set    |
| 6   | SMART Touch (Control room)               | 1 set    |

## 2.4 Specifications

### 2.4.1 Image Processing

| No. | Item                        | Specifications  |
|-----|-----------------------------|---|
| 1   | Image input                 | <ul style="list-style-type: none"> <li>Density resolution: 12 bits/4,096 steps</li> <li>Sampling frequency: 40.0 MHz</li> </ul> |
| 2   | Image operation performance | The 1024x1024 image processing (noise reduction filter processing, etc.) is performed at a maximum of 30 frames/s.b             |
| 3   | Display D/A converter       | 8 bits/256 steps or more  |

### 2.4.2 Image Display

| No. | Item                 | Specifications  |
|-----|----------------------|---|
| 1   | Output video signal  | Non-interlaced  |
| 2   | Display image matrix | 1280x1024   |
| 3   | Display step:        | 256 steps   |
| 4   | Maximum brightness   | 400 cd/m <sup>2</sup> or higher in examination room   |
| 5   | Image processing     | <ul style="list-style-type: none"> <li>Noise reduction processing</li> <li>Window level adjustment</li> <li>Edge enhancement processing</li> <li>Negative/Positive inversion</li> <li>Gamma correction</li> <li>Image magnification</li> <li>Re-masking</li> <li>Re-registration</li> <li>Virtual collimation</li> <li>Landmarking</li> <li>Peak hold</li> <li>Contour enhancement</li> </ul> |
| 6   | Character display    | <ul style="list-style-type: none"> <li>Fixed format (patient information etc.)</li> <li>Free format (annotation etc.)</li> </ul>  |
| 7   | Image analysis       | <ul style="list-style-type: none"> <li>QCA (Quantitative Coronary Analysis)</li> <li>LV (Left Ventricular Analysis)</li> </ul>  |

### 2.4.3 Image Recording

---

| No. | Item          | Specifications   |
|-----|---------------|--|
| 1   | Magnetic disk | <ul style="list-style-type: none"><li>• 1024 x 1024 pixels (12-bit): 100,000 frames</li></ul>  |
| 2   | Media         | CD-R (650 MB) <ul style="list-style-type: none"><li>• 512 x 512 pixels (8-bit): 4,800 frames or less per disk</li><li>• 1024 x 1024 pixel (12-bit): 600 frames or less per disk</li></ul> DVD-R (4.7 GB) <ul style="list-style-type: none"><li>• 512 x 512 pixel (8-bit): 30,000 frames or less per disk</li><li>• 1024 x 1024 pixel (12-bit): 4,000 frames or less per disk</li></ul> |

### 2.4.4 Network

---

| No. | Item                           |
|-----|--------------------------------|
| 1   | DICOM image/RDSR storage       |
| 2   | DICOM image/RDSR receive       |
| 3   | DICOM image storage commitment |
| 4   | DICOM image print              |
| 5   | DICOM modality worklist        |



## 2.5 Function

### 2.5.1 Image Acquisition

2

| No. | Item                   | Functions  |
|-----|------------------------|--|
| 1   | DA Mode                | Serial radiography at up to 30 frames/second is possible (for 12-bit image acquisition). Acquired images can be reviewed immediately. For example, it is easy to observe how contrast agent flows in a patient's body. Acquisition with Precession and Pendulum C-arm motion is also supported.    |
| 2   | DSA Mode               | DSA (Digital Subtraction Angiography) acquisition at up to 15 frames/second (in case of SFD-1612AF/SFD-1212AF/SFD-0808/SFD-0808AF)/12 frames (in case of SFD-1212) is possible (for 12-bit image acquisition). Acquired images can be reviewed immediately.  |
| 3   | RSM-DSA Mode           | RSM-DSA (Realtime Smoothed Mask DSA) filters enable enhanced imaging. As an alternative to DSA imaging, the RSM-DSA filter provides continuous automatic mask creation, enabling accurate enhanced imaging even with movement of the catheterization table or limb.                                |
| 4   | ROT-DA Mode            | Perform DA radiography while rotating the C-arm.   |
| 5   | RO-DSA Mode            | Perform DSA radiography while rotating the C-arm.  |
| 6   | ROT-RSM-DSA Mode       | Perform RSM-DSA radiography while rotating the C-arm.  |
| 7   | SPOT Mode              | Perform radiography of a single image.   |
| 8   | Stage Acquisition Mode | <p>For every setting frame, perform DSA acquisition while reducing the acquisition rate.</p> <p>Following acquisition rate can be used for stage acquisition mode.</p> <p>Single Radiography: 6 fps, 4 fps, 3 fps, 2 fps, 1 fps</p> <p>Bi-plane Radiography: 6 fps, 4 fps, 3 fps, 2 fps, 1 fps</p> |

### 2.5.2 Fluoroscopy

| No. | Item               | Functions   |
|-----|--------------------|---|
| 1   | Fluoroscopy Record | Images can be saved during or after fluoroscopy.  |
| 2   | FluoroMAP          | <p>FluoroMAP assists in the positioning of a catheter in a complicated blood-vessel route. It highlights in white, the path of contrast agent injection. The path is retained onscreen like a road map during subsequent Fluoro acquisitions to make it easy to see the moving catheter with little or no additional contrast agent.</p> <p>Subtract the background image such as bone, if necessary.</p> <p>In case of LIDM, a fluoroscopy image during radiation and vascular MAP image can be displayed at the same time.</p>  |
| 3   | DSA-MAP            | <p>DSA-MAP assists in the positioning of a catheter in a complicated blood-vessel route. Acquired blood-vessel image and fluoroscopy image can be displayed as a superimposed image. As it is not necessary to acquire the blood-vessel image again, radiation dosage and the volume of contrast agent can be reduced even more than FluoroMAP.</p> <p>Cannot use a blood-vessel image which is saved as a still image.</p> <p>Subtract the background image such as bone, if necessary.</p> <p>In case of LIDM, a fluoroscopy image during radiation and vascular MAP image can be displayed at the same time.</p> |
| 4   | TraceMAP           | <p>Fluoroscopy image and vessel contour image generated with DSA mode can be displayed as a superimposed. Interior of a blood-vessel is displayed in a transparent mode, so this is available when treat a large vessel.</p> <p>Cannot use a vessel contour image which is saved as a still image.</p>  |
| 5   | BlankMAP           | Subtract the background image such as bone, and display the fluoroscopy image. It is useful for checking device and medical agent added after BlankMAP started.   |

### 2.5.3 Image Processing

| No. | Item                                      | Functions   |
|-----|---|---|
| 1   | Serial Animation Display (Cyclic Display) | Serial animation up to 30 frames/second is available.   |
| 2   | Real-time Edge Enhancement                | Sharp images are acquired by emphasizing the subject edge with spatial frequency emphasis. It is selectable as either 5x5 Convolution processing or unsharp processing, as a technique of the spatial frequency emphasis. |
| 3   | Negative/Positive Inversion Display       | Images can be displayed as acquired (positive) or the inverse (negative).   |
| 4   | Gamma Correction                          | The display gamma can be selected.  |
| 5   | Animation Zoom                            | Still images and loops can be zoomed to up to 2.5 times of original size.   |
| 6   | Auto Window Control                       | Images are displayed with automatically-controlled stabilized contrast.   |
| 7   | Noise Reduction                           | Reduces image noise.  |
| 8   | Re-masking of DSA Images                  | A new sub image is created by re-selecting the mask image arbitrarily.  |
| 9   | Re-registration of DSA Images             | The mis-registration is reduced by moving the mas image up, down, right and left, and by executing subtraction at the position with the live image.   |
| 10  | Virtual Collimation                       | Displays the position of collimator on the window when moving the collimator.   |
| 11  | Landmarking                               | Modifies the Mask weight on DSA image.  |
| 12  | Peak Hold                                 | Extract and display the position of white and black peak in the image.  |
| 13  | Contour Enhancement                       | Automatically enhance the blood-vessel contour in the image.  |
| 14  | Flex-APS (Option)                         | Adjusts the position of each region of mask image on DSA image.   |

### 2.5.4 Image Analysis

| No. | Item                                 | Functions  |
|-----|--------------------------------------|--|
| 1   | QCA (Quantitative Coronary Analysis) | The level of the blood-vessel stensis is quantitatively analyzed. The vessel wall is recognized automatically when the stenosis part is specified, and the stenosis rate is displayed. |
| 2   | LV (Left Ventricular Analysis)       | The function of the left ventricle (LV) is quantitatively analyzed. When inner wall of LV is specified, LV lumen capacity and ejection fraction, etc. are calculated.                  |

### 2.5.5 Network

---

| No. | Item                           | Functions  |
|-----|--------------------------------|--|
| 1   | DICOM Image/RDSR Storage       | Image and RDSR are automatically transferred to the server via DICOM network.  |
| 2   | DICOM Image/RDSR Receive       | The server can be queried and then images and RDSR can be received via DICOM network.  |
| 3   | DICOM Image Storage Commitment | Images are transferred to the image server via DICOM network, with full Storage Commitment support. Once the Storage Commitment feature verifies that the images were successfully received by the image server, the images can be safely deleted from the local system. |
| 4   | DICOM Image Print              | Images can be printed on DICOM printers.   |
| 5   | DICOM Modality Worklist        | Study information, such as patient name, can be acquired from a worklist server. This speeds up patient information entry when creating a new study.   |

### 2.5.6 Image Recording

---

| No. | Item  | Functions   |
|-----|-------|---|
| 1   | CD-R  | The study including still images and loops can be saved on CD-R or DVD-R. |
| 2   | DVD-R |   |

### 2.5.7 C-arm Reposition

---

| No. | Item                               | Functions   |
|-----|------------------------------------|---|
| 1   | Filtering Images by C-arm Position | Display the image that matches the C-arm angle on the image selector.                                       |
| 2   | C-arm Repositioning                | The C-arm can be rotated to the same angle as the current image that is displayed on the Reference monitor. |

### 2.5.8 Table Reposition

---

| No. | Item  | Functions  |
|-----|---|--|
| 1   | Table Repositioning<br>(Combined with KS-100) | The table can be moved to the same position as the current image that is displayed on the Reference monitor. |

## 2.6 Special Information

Special information of classification, labeling and environmental conditions are described in this section.

2

### 2.6.1 Classification

#### ■ The classification of protection against electric shock:

- Class I equipment



It means that CLASS I EQUIPMENT, that is, electrical equipment in which protection against electric shock does not rely on BASIC INSULATION only, but which includes an additional safety precaution in that means are provided for ACCESSIBLE PARTS of metal or internal parts of metal to be PROTECTIVELY EARTHED. (International standard IEC 60601-1:2005)

#### ■ The degree of protection against electric shock:

- Equipment not including applied Part.

#### ■ The classification of EMC (Electro-Magnetic Compatibility):

- Group 1 and Class A



**This system belongs to Group 1 and Class A equipment in accordance with EN60601-1-2:2007.**

Class A equipment is equipment suitable for use in all establishments other than domestic and those directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.

Group 1 contains all ISM equipment in which there is intentionally generated or used conductively coupled RF energy that is necessary for the internal functioning of the equipment itself.

### ■ The degree of protection against effects of water:

- Ordinary equipment:



**This equipment is not protected against immersion in liquid.**

Never use this equipment in a place in which immersion in liquid may occur. Never spill liquid on the surface or inside the equipment. Otherwise, electrical shock may occur. When liquid is spilled, contact our service office or a Shimadzu service representative.

### ■ The degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide:

- This equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.





**Risk of explosion:**

There is a risk of explosion if the equipment is used near flammable anesthetics.

## 2.6.2 Labeling


### Warning/Caution Labels

Warning/Caution labels are attached to the side of the main system cabinet like this:



## 2 System Overview

Warning label is attached to the System Display like this:



Caution and unit labels are attached to the back and side of water cooling unit like this:

(1) HEC002-A5B-X101



(2) HECR002-A5






## Storage and Transport Labels

The label for storage and transport condition is attached on the package.




## Others


Labels for Digital Acquisition Cabinet are located as follows:



(1) Rating Label




OR



(2) cTUVus



### (3) FDA Certification



The rating label for Digital Acquisition Cabinet Lateral (option) is located as follows:




|   |                   |    |       |
|---|-------------------|----|-------|
| MODEL   | ACQ3 CABINET(L)   |    |       |
| SERIAL NO.  |                   |    |       |
| SUPPLY MAINS  |                   |    |       |
| V   | ~ 200/220/230/240 | Hz | 50/60 |
| INPUT POWER   | kVA               | 3  |       |
| MANUFACTURED:   |                   |    | KYOTO |
| PROTECTION CLASS  | I                 |    |       |
| 555-11395   |                   |    |       |
| 株式会社島津製作所 604-8511<br>京都市中京区西ノ京桑原町1<br><b>SHIMADZU CORPORATION</b><br>1, Nishinokyo-Kuwaeracho, Nakagyo-ku, Kyoto 604-8511, Japan |                   |    |       |

The rating label for Transformer Box (option) is located as follows:




|   |                   |      |       |
|---|-------------------|------|-------|
| MODEL   | TRANS BOX ASSY    |      |       |
| SERIAL NO.  | 018S134001        |      |       |
| SUPPLY MAINS  |                   |      |       |
| V   | ~ 200/220/230/240 | Hz   | 50/60 |
| INPUT POWER   | kVA               | 3    |       |
| MANUFACTURED:   | APRIL             | 2006 | KYOTO |
| PROTECTION CLASS  | I                 |      |       |
| 503-45621   |                   |      |       |
| 株式会社島津製作所 604-8511<br>京都市中京区西ノ京桑原町1<br><b>SHIMADZU CORPORATION</b><br>1, Nishinokyo-Kuwaeracho, Nakagyo-ku, Kyoto 604-8511, Japan |                   |      |       |

The rating label for Large Monitor Cabinet (option) is located as follows:





2

The rating label of Large Monitor Cabinet 2 (option) is located as follows:




The rating label for DVI OUT BOX (option) is located as follows:




## 2 System Overview

The rating label for SMART Touch is located as follows:



















The rating label for Wireless Mouse is located as follows:



### 2.6.3 Symbols

Symbols used on this system are shown and described as follow:

| No. | Symbol  | Location  | Meaning  |
|-----|---|---|--|
| 1   |    | On the back of the Digital processing system cabinet.   | Risk of fire.  |
| 2   |    | On warning/caution labels.  | Attention: refer to Operation Manual and attachment  |
| 3   |    | On System display   | Turns on the equipment power.                        |
| 4   |    | On System display.  | Turns off the equipment power.                       |
| 5   |    | On Digital processing system cabinet.   | Turns on or off the equipment power.                 |
| 6   |    | On the breaker switch on the back of the Digital processing system cabinet and the transformer box. | Turns on the equipment power.                        |
| 7   |   | On the breaker switch on the back of the Digital processing system cabinet and the transformer box. | Turns off the equipment power.                       |
| 8   |  | Caution labels.   | Alternating current.                                 |
| 9   |  | Protective earth of power cable inside the system.  | Protective earth.                                    |
| 10  |  | In name plate on covers.  | Refer to Operation Manual.                           |
| 11  |  | In name plate on covers.  | Date of manufacture.                                 |
| 12  |  | In name plate on covers.  | Manufacturer.  |
| 13  |  | In name plate on covers.  | Serial number.                                       |
| 14  |  | On the foot stand of the System display.  | Ionizing radiation.                                  |
| 15  |  | Integrated acquisition cabinet.   | Authorized Representative in the European Community. |
| 15  |  | On the USB cable for wireless mouse.  | RF transmitter.                                      |

## 2.6.4 Environmental Conditions

### ■ Setup and Operating

| No. | Item                    | Environmental Conditions   |
|-----|-------------------------|--|
| 1   | Power Source            | Digital Acquisition Cabinet:<br>Single phase AC200/220/230/240 V, 3 kVA, 50/60 Hz<br>Transformer Box *:<br>Single phase AC 200/220/230/240 V, 3 kVA, 50/60 Hz<br>A power source that is different than that of C-arm and X-ray High Voltage Control unit is required.<br>(* )Used when adding an optional computer such as a 3D Workstation. |
| 2   | Temperature/Humidity    | <ul style="list-style-type: none"> <li>Examination Room<br/>from 10 °C to 35 °C<br/>from 15 % to 75 % (non condensing)</li> <li>Control Room<br/>from 10 °C to 35 °C<br/>from 15 % to 75 % (non condensing)</li> <li>Equipment Room<br/>from 10 °C to 30 °C<br/>from 15 % to 75 % (non condensing)</li> </ul>                                |
| 3   | Grounding               | Grounding resistance: 100 Ohm max.   |
| 4   | Generated Heat Quantity | <ul style="list-style-type: none"> <li>Equipment Room<br/>Digital acquisition cabinet: 1560 Kcal/h<br/>FPD cooling system: 1000 Kcal/h</li> <li>Control Room<br/>Digital reference cabinet: 470 Kcal/h<br/>LCD Monitor 48 Kcal/h</li> </ul>  |



- When isolating the equipment from the power supply, open the circuit breaker or knife switch of the switchboard that the equipment is connected to.
- Prepare some locking mechanism to keep an OFF position of the circuit breaker and knife switch.

### ■ Storage and Transportation

| No. | Item        | Environmental Conditions |
|-----|-------------|--------------------------|
| 1   | Temperature | from -10 °C to 60 °C     |
| 2   | Humidity    | from 10 °C to 95 °C      |
| 3   | Pressure    | from 700 hPa to 1060 hPa |



## ■ Conformity of Equipment

Refer to the Operation Manual which system to be combined.

### 2.6.5 EMC (Electro Magnetic Compatibility) Information

#### ■ Operating Condition

This system belongs to Group 1 and Class A equipment in accordance to IEC60601-1-2:2007.

The system uses radio-frequency energy only for its internal function and is not intended to deliver energy to the patient. But little leakage of radio-frequency energy does harm to high-sensitive equipment.

The system main power line in the clinical site should be connected to the domestic power source which are separated from the public main network.

DAR-9500f needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.

Portable and mobile RF communications equipment can affect the DAR-9500f.

### 2.6.6 Wireless Communications

#### ■ Wireless Mouse

| Item | Frequency Range | Modulation | Output Watts |
|------|-----------------|------------|--------------|
| -    | 2400-2483 MHz   | GFSK       | Max: 1.5 mW  |

### Cable List

- In case of ACQ3

| Digital Angiography System DAR-9500f |                                   |              |                  |          |
|--------------------------------------|-----------------------------------|--------------|------------------|----------|
| No.                                  | Cable Type                        | Manufacturer | Cable Length (m) | Shielded |
| 1                                    | AC Cable2 CCC                     | SHIMADZU     | 10               | No       |
| 2                                    | Earth Cable                       | SHIMADZU     | 10               | No       |
| 3                                    | AC DAR-TR Cable                   | SHIMADZU     | 20               | No       |
| 4                                    | Earth Cable                       | SHIMADZU     | 20               | No       |
| 5                                    | LAN Cable                         | SHIMADZU     | 20               | No       |
| 6                                    | Optical Power Cable               | SHIMADZU     | 40               | Yes      |
| 7                                    | Power Cable attached with Chiller | SHIMADZU     | 10               | No       |
| 8                                    | Cooler Warning Cable              | SHIMADZU     | 10               | No       |
| 9                                    | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 10                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 11                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 12                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 13                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 14                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 15                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 16                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 17                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 18                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 19                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 20                                   | LAN DVI Extension Cable           | SHIMADZU     | 30               | Yes      |
| 21                                   | DVI PS Cable                      | SHIMADZU     | 30               | Yes      |
| 22                                   | DVI PS Cable                      | SHIMADZU     | 30               | Yes      |
| 23                                   | DVI PS Cable                      | SHIMADZU     | 30               | Yes      |
| 24                                   | AC Power Cable                    | SHIMADZU     | 2                | No       |
| 25                                   | AC Power Cable                    | SHIMADZU     | 2                | No       |
| 26                                   | AC Power Cable                    | SHIMADZU     | 2                | No       |
| 27                                   | AC Power Cable                    | SHIMADZU     | 2                | No       |
| 28                                   | DVI Cable                         | SHIMADZU     | 2                | Yes      |
| 29                                   | DVI Cable                         | SHIMADZU     | 2                | Yes      |



| Digital Angiography System DAR-9500f |                       |              |                  |          |
|--------------------------------------|-----------------------|--------------|------------------|----------|
| No.                                  | Cable Type            | Manufacturer | Cable Length (m) | Shielded |
| 30                                   | DVI Cable             | SHIMADZU     | 2                | Yes      |
| 31                                   | DVI Cable             | SHIMADZU     | 2                | Yes      |
| 32                                   | AC Power Cable        | SHIMADZU     | 2                | No       |
| 33                                   | AC Power Cable        | SHIMADZU     | 2                | No       |
| 36                                   | DVI Cable             | SHIMADZU     | 5                | Yes      |
| 37                                   | DVI Cable             | SHIMADZU     | 5                | Yes      |
| 40                                   | IVR NEO Cable ASSY    | SHIMADZU     | 20               | Yes      |
| 42                                   | Power Cable Joint Box | SHIMADZU     | 5                | No       |
| 43                                   | Earth Cable Joint Box | SHIMADZU     | 5                | No       |
| 44                                   | IVR J11 Cable         | SHIMADZU     | 20               | Yes      |
| 45                                   | JOINT BOX Cable       | SHIMADZU     | 20               | Yes      |
| 46                                   | VGA Cable             | SHIMADZU     | 30               | Yes      |
| 47                                   | VGA Cable             | SHIMADZU     | 30               | Yes      |
| 48                                   | AUDIO SD Cable        | SHIMADZU     | 30               | Yes      |
| 49                                   | Audio ID Cable        | SHIMADZU     | 30               | Yes      |
| 50                                   | AC Power Cable        | SHIMADZU     | 3                | No       |
| 51                                   | Earth Cable           | SHIMADZU     | 3                | No       |
| 52                                   | SYS-DISP Cable        | SHIMADZU     | 3                | Yes      |
| 53                                   | Touch Panel I/F Cable | SHIMADZU     | 3                | Yes      |
| 54                                   | AC Power Cable        | SHIMADZU     | 2                | No       |
| 55                                   | IVR J10 Cable         | SHIMADZU     | 3                | Yes      |
| 56                                   | J1F Cable             | SHIMADZU     | 20               | Yes      |
| 57                                   | RSYS-XGENT-F Cable    | SHIMADZU     | 10               | Yes      |
| 58                                   | DAR-ARCNET Cable      | SHIMADZU     | 10               | Yes      |
| 59                                   | IBS Cable             | SHIMADZU     | 20               | Yes      |
| 60                                   | Serial SDPC           | SHIMADZU     | 15               | Yes      |
| 61                                   | SA60 Serial Cable     | SHIMADZU     | 10               | Yes      |
| 62                                   | RS Cable              | SHIMADZU     | 15               | Yes      |
| 63                                   | DOSE Meter Cable      | SHIMADZU     | 30               | Yes      |
| 64                                   | Modular Cable         | SHIMADZU     | 20               | No       |
| 66                                   | LAN Cable (cross)     | SHIMADZU     | 15               | No       |

## 2 System Overview

| Digital Angiography System DAR-9500f |                        |              |                  |          |
|--------------------------------------|------------------------|--------------|------------------|----------|
| No.                                  | Cable Type             | Manufacturer | Cable Length (m) | Shielded |
| 67                                   | AC Power Cable         | SHIMADZU     | 2.5              | No       |
| 68                                   | Cont. Microphone Cable | SHIMADZU     | 3                | Yes      |
| 69                                   | Cont. Microphone Cable | SHIMADZU     | 5                | Yes      |
| 70                                   | Cont. Speaker Cable    | SHIMADZU     | 10               | No       |
| 71                                   | Exam. Microphone Cable | SHIMADZU     | 30               | Yes      |
| 72                                   | Exam. Speaker Cable    | SHIMADZU     | 30               | Yes      |
| 73                                   | BNC Cable              | SHIMADZU     | 30               | Yes      |
| 74                                   | AC Cable2 CCC          | SHIMADZU     | 10               | No       |
| 75                                   | Earth Cable            | SHIMADZU     | 10               | No       |
| 76                                   | AC Cable               | SHIMADZU     | 20               | No       |
| 77                                   | Earth Cable            | SHIMADZU     | 20               | No       |
| 78                                   | AC Cable               | SHIMADZU     | 20               | No       |
| 79                                   | Earth Cable            | SHIMADZU     | 20               | No       |
| 80                                   | LAN Cable              | SHIMADZU     | 30               | No       |
| 81                                   | LAN Cable              | SHIMADZU     | 30               | No       |
| 82                                   | LAN Cable              | SHIMADZU     | 30               | No       |
| 83                                   | LAN Cable              | SHIMADZU     | 30               | No       |
| 84                                   | AC Power Cable         | SHIMADZU     | 2                | No       |
| 85                                   | AC Power Cable         | SHIMADZU     | 2                | No       |
| 86                                   | AC Power Cable         | SHIMADZU     | 2                | No       |
| 87                                   | AC Power Cable         | SHIMADZU     | 2                | No       |
| 88                                   | AC Power Cable         | SHIMADZU     | 2                | No       |
| 89                                   | AC Power Cable         | SHIMADZU     | 2                | No       |
| 90                                   | AC Power Cable         | SHIMADZU     | 2                | No       |
| 91                                   | AC Power Cable         | SHIMADZU     | 2                | No       |
| 92                                   | DVI Cable              | SHIMADZU     | 5                | Yes      |
| 94                                   | DVI Cable              | SHIMADZU     | 5                | Yes      |
| 95                                   | DVI Cable              | SHIMADZU     | 5                | Yes      |
| 96                                   | DVI Cable              | SHIMADZU     | 5                | Yes      |
| 97                                   | DVI Cable              | SHIMADZU     | 5                | Yes      |
| 98                                   | DVI Cable              | SHIMADZU     | 5                | Yes      |

| Digital Angiography System DAR-9500f |                           |              |                  |          |
|--------------------------------------|---------------------------|--------------|------------------|----------|
| No.                                  | Cable Type                | Manufacturer | Cable Length (m) | Shielded |
| 99                                   | DVI Cable                 | SHIMADZU     | 5                | Yes      |
| 100                                  | DVI Cable                 | SHIMADZU     | 5                | Yes      |
| 101                                  | DVI Cable                 | SHIMADZU     | 5                | Yes      |
| 102                                  | DVI Cable                 | SHIMADZU     | 5                | Yes      |
| 103                                  | DVI Cable                 | SHIMADZU     | 5                | Yes      |
| 104                                  | VGA Cable                 | SHIMADZU     | 30               | Yes      |
| 105                                  | VGA Cable                 | SHIMADZU     | 30               | Yes      |
| 106                                  | VGA Cable                 | SHIMADZU     | 30               | Yes      |
| 107                                  | VGA Cable                 | SHIMADZU     | 30               | Yes      |
| 108                                  | KVM DVI Cable             | SHIMADZU     | 1.8              | Yes      |
| 109                                  | KVM DVI Cable             | SHIMADZU     | 1.8              | Yes      |
| 110                                  | DVI Cable                 | SHIMADZU     | 5                | Yes      |
| 111                                  | AC power Cable            | SHIMADZU     | 1.9              | No       |
| 112                                  | KVM VGA Cable             | SHIMADZU     | 1.2              | Yes      |
| 113                                  | KVM VGA Cable             | SHIMADZU     | 1.2              | Yes      |
| 114                                  | Mouse EXT Cable           | SHIMADZU     | 20               | Yes      |
| 115                                  | Earth Cable               | SHIMADZU     | 20               | No       |
| 116                                  | AC Power Cable            | SHIMADZU     | 2                | No       |
| 117                                  | AC Power Cable            | SHIMADZU     | 2                | No       |
| 118                                  | DVI Cable                 | SHIMADZU     | 5                | Yes      |
| 119                                  | DVI Cable                 | SHIMADZU     | 5                | Yes      |
| 120                                  | DVI Cable                 | SHIMADZU     | 5                | Yes      |
| 121                                  | LAN Cable                 | SHIMADZU     | 30               | No       |
| 122                                  | G40 Serial Cable          | SHIMADZU     | 20               | Yes      |
| 123                                  | LAN Cable                 | SHIMADZU     | 30               | No       |
| 125                                  | DVI LAN CABLE             | SHIMADZU     | 36               | Yes      |
| 126                                  | DVI LAN CABLE             | SHIMADZU     | 36               | Yes      |
| 127                                  | DVI LAN CABLE             | SHIMADZU     | 36               | Yes      |
| 128                                  | DVI LAN CABLE             | SHIMADZU     | 36               | Yes      |
| 129                                  | AC Power Cable            | SHIMADZU     | 2                | No       |
| 130                                  | Power Cable Large Monitor | SHIMADZU     | 30               | Yes      |

## 2 System Overview

| Digital Angiography System DAR-9500f |                         |              |                  |          |
|--------------------------------------|-------------------------|--------------|------------------|----------|
| No.                                  | Cable Type              | Manufacturer | Cable Length (m) | Shielded |
| 131                                  | Earth Cable             | SHIMADZU     | 30               | No       |
| 132                                  | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 133                                  | DVI Extension Cable     | SHIMADZU     | 36               | Yes      |
| 134                                  | Power Cable             | SHIMADZU     | 30               | No       |
| 135                                  | Mouse                   | SHIMADZU     | 2                | No       |
| 136                                  | Keyboard                | SHIMADZU     | 2                | No       |
| 137                                  | Mouse                   | SHIMADZU     | 2                | No       |
| 138                                  | Keyboard                | SHIMADZU     | 2                | No       |
| 139                                  | Mouse                   | SHIMADZU     | 2                | No       |
| 140                                  | Keyboard                | SHIMADZU     | 2                | No       |
| 141                                  | Mouse                   | SHIMADZU     | 2                | No       |
| 142                                  | Keyboard                | SHIMADZU     | 2                | No       |
| 143                                  | Mouse                   | SHIMADZU     | 2                | No       |
| 144                                  | Keyboard                | SHIMADZU     | 2                | No       |
| 145                                  | Mouse                   | SHIMADZU     | 2                | No       |
| 146                                  | Keyboard                | SHIMADZU     | 2                | No       |
| 147                                  | Power Cable             | SHIMADZU     | 30               | No       |
| 148                                  | Power Cable             | SHIMADZU     | 30               | No       |
| 149                                  | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 150                                  | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 151                                  | AC Cable2 CCC           | SHIMADZU     | 10               | No       |
| 152                                  | Earth Cable             | SHIMADZU     | 10               | No       |
| 153                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 154                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 155                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 156                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 157                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 158                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 159                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 160                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 161                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |

| Digital Angiography System DAR-9500f |                         |              |                  |          |
|--------------------------------------|-------------------------|--------------|------------------|----------|
| No.                                  | Cable Type              | Manufacturer | Cable Length (m) | Shielded |
| 162                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 163                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 164                                  | LAN DVI Extension Cable | SHIMADZU     | 30               | Yes      |
| 165                                  | DVI PS Cable            | SHIMADZU     | 30               | Yes      |
| 166                                  | DVI PS Cable            | SHIMADZU     | 30               | Yes      |
| 167                                  | DVI PS Cable            | SHIMADZU     | 30               | Yes      |
| 168                                  | AC Power Cable          | SHIMADZU     | 2                | No       |
| 169                                  | AC Power Cable          | SHIMADZU     | 2                | No       |
| 170                                  | AC Power Cable          | SHIMADZU     | 2                | No       |
| 171                                  | AC Power Cable          | SHIMADZU     | 2                | No       |
| 172                                  | DVI Cable               | SHIMADZU     | 2                | Yes      |
| 173                                  | DVI Cable               | SHIMADZU     | 2                | Yes      |
| 174                                  | DVI Cable               | SHIMADZU     | 2                | Yes      |
| 175                                  | DVI Cable               | SHIMADZU     | 2                | Yes      |
| 176                                  | AC Power Cable          | SHIMADZU     | 2                | No       |
| 177                                  | AC Power Cable          | SHIMADZU     | 2                | No       |
| 178                                  | DVI Cable               | SHIMADZU     | 2                | Yes      |
| 179                                  | DVI Cable               | SHIMADZU     | 2                | Yes      |
| 180                                  | SOLIOS CABLE            | SHIMADZU     | 10               | Yes      |
| 181                                  | RSYS-CONV CABLE         | SHIMADZU     | 10               | Yes      |
| 182                                  | LAN Cable               | SHIMADZU     | 10               | Yes      |
| 183                                  | AC Power Cable          | SHIMADZU     | 10               | No       |
| 184                                  | RS232C CABLE            | SHIMADZU     | 10               | Yes      |
| 185                                  | LAN DVI Cable           | SHIMADZU     | 30               | Yes      |
| 186                                  | Relay Cable             | SHIMADZU     | 10               | Yes      |
| 187                                  | Optical Power Cable     | SHIMADZU     | 40               | Yes      |
| 188                                  | Power Cable for Chiller | SHIMADZU     | 10               | No       |
| 189                                  | Cooler Warning Cable    | SHIMADZU     | 10               | No       |
| 190                                  | RSYS-XGEN-L Cable       | SHIMADZU     | 20               | Yes      |
| 191                                  | IBS Cable               | SHIMADZU     | 20               | Yes      |
| 192                                  | Serial SDPC             | SHIMADZU     | 20               | Yes      |

## 2 System Overview

| Digital Angiography System DAR-9500f |                     |              |                  |          |
|--------------------------------------|---------------------|--------------|------------------|----------|
| No.                                  | Cable Type          | Manufacturer | Cable Length (m) | Shielded |
| 193                                  | G40 Serial Cable    | SHIMADZU     | 20               | Yes      |
| 194                                  | J1L Cable           | SHIMADZU     | 20               | Yes      |
| 195                                  | AC Power Cable      | SHIMADZU     | 10               | No       |
| 196                                  | LAN Cable           | SHIMADZU     | 30               | Yes      |
| 197                                  | DOSE Meter Cable    | SHIMADZU     | 30               | Yes      |
| 198                                  | AC Cable            | SHIMADZU     | 20               | No       |
| 199                                  | Earth Cable         | SHIMADZU     | 20               | No       |
| 200                                  | LAN Cable           | SHIMADZU     | 30               | Yes      |
| 201                                  | DVI Extension Cable | SHIMADZU     | 36               | Yes      |
| 202                                  | AC Power Cable      | SHIMADZU     | 2.5              | No       |
| 203                                  | AC Power Cable      | SHIMADZU     | 3                | No       |
| 204                                  | USB Extension Cable | SHIMADZU     | 30               | Yes      |
| 205                                  | Keyboard            | SHIMADZU     | 2                | No       |
| 206                                  | Mouse               | SHIMADZU     | 2                | No       |
| 207                                  | BNC Cable           | SHIMADZU     | 30               | Yes      |
| 208                                  | UPS Serial Cable    | SHIMADZU     | 20               | Yes      |
| 209                                  | MH-400 Sync Cable   | SHIMADZU     | 20               | Yes      |
| 219                                  | DVI LAN CABLE       | SHIMADZU     | 36               | Yes      |
| 220                                  | AC POWER CABLE      | SHIMADZU     | 2                | No       |
| 221                                  | LAN CABLE           | SHIMADZU     | 30               | No       |
| 222                                  | LAN CABLE           | SHIMADZU     | 30               | No       |
| 223                                  | LAN CABLE           | SHIMADZU     | 30               | No       |
| 224                                  | LAN CABLE           | SHIMADZU     | 30               | No       |
| 225                                  | LAN CABLE           | SHIMADZU     | 30               | No       |
| 226                                  | LAN CABLE           | SHIMADZU     | 30               | No       |
| 227                                  | LAN CABLE           | SHIMADZU     | 30               | No       |
| 228                                  | LAN CABLE           | SHIMADZU     | 30               | No       |
| 230                                  | Hand Switch Cable   | SHIMADZU     | 0.2              | No       |
| 231                                  | AC Power Cable      | SHIMADZU     | 3                | No       |
| 232                                  | Earth Cable         | SHIMADZU     | 3                | No       |
| 233                                  | USB Cable           | SHIMADZU     | 0.9              | Yes      |

| Digital Angiography System DAR-9500f |                         |              |                  |          |
|--------------------------------------|-------------------------|--------------|------------------|----------|
| No.                                  | Cable Type              | Manufacturer | Cable Length (m) | Shielded |
| 234                                  | AC Power Cable          | SHIMADZU     | 30               | No       |
| 235                                  | Power Cable Touch Panel | SHIMADZU     | 4                | Yes      |
| 236                                  | Earth Cable             | SHIMADZU     | 4                | No       |
| 237                                  | USB Cable               | SHIMADZU     | 5                | Yes      |
| 238                                  | Power Cable Touch Panel | SHIMADZU     | 2                | Yes      |
| 239                                  | Earth Cable             | SHIMADZU     | 2                | No       |
| 240                                  | USB Cable               | SHIMADZU     | 0.9              | Yes      |
| 241                                  | LAN CABLE               | SHIMADZU     | 30               | Yes      |
| 242                                  | LAN CABLE               | SHIMADZU     | 30               | Yes      |
| 243                                  | LAN CABLE               | SHIMADZU     | 30               | Yes      |
| 244                                  | USB EXT Cable           | SHIMADZU     | 30               | Yes      |
| 245                                  | Earth Cable             | SHIMADZU     | 0.9              | No       |
| 246                                  | AC Power Cable          | SHIMADZU     | 30               | No       |
| 247                                  | LAN CABLE               | SHIMADZU     | 30               | No       |
| 248                                  | LAN CABLE               | SHIMADZU     | 30               | No       |
| 249                                  | AC Power Cable          | SHIMADZU     | 2                | No       |
| 250                                  | AC Power Cable          | SHIMADZU     | 2                | No       |
| 251                                  | DC Power Cable          | SHIMADZU     | 0.5              | No       |
| 252                                  | LAN CABLE               | SHIMADZU     | 30               | Yes      |
| 253                                  | DVI CABLE               | SHIMADZU     | 5                | Yes      |
| 254                                  | DVI CABLE               | SHIMADZU     | 5                | Yes      |
| 255                                  | Keyboard                | SHIMADZU     | 4                | No       |
| 256                                  | Mouse                   | SHIMADZU     | 4                | No       |
| 257                                  | USB Cable               | SHIMADZU     | 5                | No       |
| 258                                  | USB Cable               | SHIMADZU     | 5                | No       |
| 259                                  | USB Cable               | SHIMADZU     | 5                | No       |
| 260                                  | LAN CABLE               | SHIMADZU     | 30               | Yes      |

## 2 System Overview

- In case of ACQ4

| Digital Angiography System DAR-9500f |                         |              |                  |          |
|--------------------------------------|-------------------------|--------------|------------------|----------|
| No.                                  | Cable Type              | Manufacturer | Cable Length (m) | Shielded |
| 1                                    | AC Cable2 CCC           | SHIMADZU     | 10               | No       |
| 2                                    | Earth Cable             | SHIMADZU     | 10               | No       |
| 3                                    | JOINT BOX CABLE         | SHIMADZU     | 25               | Yes      |
| 4                                    | DOSE Meter Cable        | SHIMADZU     | 35               | Yes      |
| 5                                    | Optical Power Cable     | SHIMADZU     | 40               | Yes      |
| 6                                    | Power Cable for Chiller | SHIMADZU     | 10               | No       |
| 7                                    | Cooler Warning Cable    | SHIMADZU     | 10               | No       |
| 10                                   | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 11                                   | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 12                                   | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 15                                   | ISYS-XGEN Cable         | SHIMADZU     | 20               | Yes      |
| 16                                   | IBS Cable               | SHIMADZU     | 20               | Yes      |
| 17                                   | ACQ IO-XGEN CABLE       | SHIMADZU     | 20               | Yes      |
| 20                                   | SYS DISP CABLE          | SHIMADZU     | 3                | Yes      |
| 21                                   | Earth Cable             | SHIMADZU     | 3                | No       |
| 22                                   | PC IF POWER CABLE       | SHIMADZU     | 3                | Yes      |
| 23                                   | Hand Switch Cable       | SHIMADZU     | 0.2              | No       |
| 24                                   | Keyboard                | SHIMADZU     | 2                | No       |
| 25                                   | Mouse                   | SHIMADZU     | 2                | No       |
| 26                                   | Earth Cable             | SHIMADZU     | 1.5              | Yes      |
| 27                                   | USB Cable               | SHIMADZU     | 1.5              | No       |
| 28                                   | DVI Cable               | SHIMADZU     | 5                | Yes      |
| 29                                   | DC Power Cable          | SHIMADZU     | 1.5              | No       |
| 30                                   | DVI Cable               | SHIMADZU     | 2                | Yes      |
| 31                                   | DVI Cable               | SHIMADZU     | 2                | Yes      |
| 32                                   | DVI Cable               | SHIMADZU     | 2                | Yes      |
| 33                                   | AC Power Cable          | SHIMADZU     | 2                | No       |
| 34                                   | AC Power Cable          | SHIMADZU     | 2                | No       |
| 35                                   | AC Power Cable          | SHIMADZU     | 2                | No       |
| 40                                   | AC Power Cable          | SHIMADZU     | 30               | No       |



| Digital Angiography System DAR-9500f |                         |              |                  |          |
|--------------------------------------|-------------------------|--------------|------------------|----------|
| No.                                  | Cable Type              | Manufacturer | Cable Length (m) | Shielded |
| 41                                   | Power Cable Touch Panel | SHIMADZU     | 4                | Yes      |
| 42                                   | Earth Cable             | SHIMADZU     | 4                | No       |
| 43                                   | USB Cable               | SHIMADZU     | 3                | Yes      |
| 44                                   | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 45                                   | AC Power Cable          | SHIMADZU     | 30               | No       |
| 46                                   | Power Cable Touch Panel | SHIMADZU     | 2                | Yes      |
| 47                                   | Earth Cable             | SHIMADZU     | 2                | No       |
| 48                                   | USB Cable               | SHIMADZU     | 0.9              | Yes      |
| 49                                   | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 50                                   | AC Power Cable          | SHIMADZU     | 3                | No       |
| 51                                   | Earth Cable             | SHIMADZU     | 3                | No       |
| 52                                   | USB Cable               | SHIMADZU     | 0.9              | Yes      |
| 53                                   | Earth Cable             | SHIMADZU     | 0.9              | No       |
| 54                                   | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 60                                   | AC Power Cable          | SHIMADZU     | 10               | No       |
| 61                                   | LAN Cable               | SHIMADZU     | 30               | Yes      |
| 62                                   | Keyboard                | SHIMADZU     | 2                | No       |
| 63                                   | Mouse                   | SHIMADZU     | 2                | No       |
| 64                                   | LAN Cable               | SHIMADZU     | 30               | No       |
| 65                                   | LAN Cable               | SHIMADZU     | 30               | No       |
| 66                                   | DVI Cable               | SHIMADZU     | 5                | Yes      |
| 67                                   | DVI Cable               | SHIMADZU     | 5                | Yest     |
| 68                                   | AC Power Cable          | SHIMADZU     | 2                | No       |
| 69                                   | AC Power Cable          | SHIMADZU     | 2                | No       |
| 70                                   | Keyboard                | SHIMADZU     | 2                | No       |
| 71                                   | Mouse                   | SHIMADZU     | 2                | No       |
| 72                                   | LAN Cable               | SHIMADZU     | 30               | No       |
| 73                                   | DVI Cable               | SHIMADZU     | 5                | Yes      |
| 74                                   | AC Power Cable          | SHIMADZU     | 2                | No       |
| 75                                   | AC Power Cable          | SHIMADZU     | 2                | No       |
| 80                                   | AC Power Cable          | SHIMADZU     | 2.5              | No       |

## 2 System Overview

| Digital Angiography System DAR-9500f |                        |              |                  |          |
|--------------------------------------|------------------------|--------------|------------------|----------|
| No.                                  | Cable Type             | Manufacturer | Cable Length (m) | Shielded |
| 81                                   | Cont. Microphone Cable | SHIMADZU     | 3                | Yes      |
| 82                                   | Cont. Microphone Cable | SHIMADZU     | 5                | Yes      |
| 83                                   | Cont. Speaker Cable    | SHIMADZU     | 10               | No       |
| 84                                   | Exam. Microphone Cable | SHIMADZU     | 30               | Yes      |
| 85                                   | Exam. Speaker Cable    | SHIMADZU     | 30               | Yes      |
| 86                                   | Cont. Microphone Cable | SHIMADZU     | 3                | Yes      |
| 87                                   | Cont. Microphone Cable | SHIMADZU     | 5                | Yes      |
| 88                                   | Cont. Speaker Cable    | SHIMADZU     | 10               | No       |
| 89                                   | Exam. Microphone Cable | SHIMADZU     | 30               | Yes      |
| 90                                   | Exam. Speaker Cable    | SHIMADZU     | 30               | Yes      |
| 91                                   | Modular Cable          | SHIMADZU     | 20               | No       |
| 92                                   | LAN Cable              | SHIMADZU     | 30               | No       |
| 95                                   | Mouse EXT Cable        | SHIMADZU     | 30               | Yes      |
| 96                                   | Earth Cable            | SHIMADZU     | 20               | No       |
| 97                                   | USB Extension Cable    | SHIMADZU     | 30               | Yes      |
| 100                                  | DVI LAN Cable          | SHIMADZU     | 36               | Yes      |
| 101                                  | DVI LAN Cable          | SHIMADZU     | 36               | Yes      |
| 102                                  | DVI LAN Cable          | SHIMADZU     | 36               | Yes      |
| 103                                  | DVI LAN Cable          | SHIMADZU     | 36               | Yes      |
| 104                                  | Power Cable LM         | SHIMADZU     | 30               | Yes      |
| 105                                  | Earth Cable            | SHIMADZU     | 30               | No       |
| 106                                  | LAN Cable              | SHIMADZU     | 30               | Yes      |
| 107                                  | DVI Extension Cable    | SHIMADZU     | 36               | Yes      |
| 110                                  | AC Cable 2 CCC         | SHIMADZU     | 10               | No       |
| 111                                  | Earth Cable            | SHIMADZU     | 10               | No       |
| 112                                  | AC Cable               | SHIMADZU     | 20               | No       |
| 113                                  | Earth Cable            | SHIMADZU     | 20               | No       |
| 114                                  | AC Cable               | SHIMADZU     | 20               | No       |
| 115                                  | Earth Cable            | SHIMADZU     | 20               | No       |
| 120                                  | AC Power Cable         | SHIMADZU     | 2                | No       |
| 121                                  | LAN Cable              | SHIMADZU     | 30               | No       |

| Digital Angiography System DAR-9500f |                   |              |                  |          |
|--------------------------------------|-------------------|--------------|------------------|----------|
| No.                                  | Cable Type        | Manufacturer | Cable Length (m) | Shielded |
| 122                                  | Mouse             | SHIMADZU     | 2                | No       |
| 123                                  | Keyboard          | SHIMADZU     | 2                | No       |
| 124                                  | DVI Cable         | SHIMADZU     | 5                | Yes      |
| 125                                  | AC Power Cable    | SHIMADZU     | 2                | No       |
| 130                                  | AC Power Cable    | SHIMADZU     | 2                | No       |
| 131                                  | LAN Cable         | SHIMADZU     | 30               | No       |
| 132                                  | Mouse             | SHIMADZU     | 2                | No       |
| 133                                  | Keyboard          | SHIMADZU     | 2                | No       |
| 134                                  | DVI Cable         | SHIMADZU     | 5                | Yes      |
| 135                                  | AC Power Cable    | SHIMADZU     | 2                | No       |
| 140                                  | AC Power Cable    | SHIMADZU     | 2                | No       |
| 141                                  | LAN Cable         | SHIMADZU     | 30               | No       |
| 142                                  | Mouse             | SHIMADZU     | 2                | No       |
| 143                                  | Keyboard          | SHIMADZU     | 2                | No       |
| 144                                  | DVI Cable         | SHIMADZU     | 5                | Yes      |
| 145                                  | AC Power Cable    | SHIMADZU     | 2                | No       |
| 150                                  | AC Power Cable    | SHIMADZU     | 2                | No       |
| 151                                  | AC Power Cable    | SHIMADZU     | 2                | No       |
| 152                                  | DC Power Cable    | SHIMADZU     | 0.5              | No       |
| 153                                  | LAN Cable         | SHIMADZU     | 30               | Yes      |
| 154                                  | DVI Cable         | SHIMADZU     | 5                | Yes      |
| 155                                  | DVI Cable         | SHIMADZU     | 5                | Yes      |
| 156                                  | Keyboard          | SHIMADZU     | 4                | No       |
| 157                                  | Mouse             | SHIMADZU     | 4                | No       |
| 158                                  | USB Cable         | SHIMADZU     | 5                | No       |
| 159                                  | USB Cable         | SHIMADZU     | 30               | No       |
| 160                                  | USB Cable         | SHIMADZU     | 5                | Yes      |
| 170                                  | LAN Cable (cross) | SHIMADZU     | 15               | No       |
| 171                                  | LAN Cable         | SHIMADZU     | 30               | No       |
| 172                                  | LAN Cable         | SHIMADZU     | 30               | No       |
| 173                                  | LAN Cable         | SHIMADZU     | 30               | No       |

## 2 System Overview

---

| Digital Angiography System DAR-9500f |                |              |                  |          |
|--------------------------------------|----------------|--------------|------------------|----------|
| No.                                  | Cable Type     | Manufacturer | Cable Length (m) | Shielded |
| 174                                  | LAN Cable      | SHIMADZU     | 30               | No       |
| 175                                  | LAN Cable      | SHIMADZU     | 30               | No       |
| 176                                  | LAN Cable      | SHIMADZU     | 30               | No       |
| 177                                  | LAN Cable      | SHIMADZU     | 30               | No       |
| 178                                  | LAN Cable      | SHIMADZU     | 30               | No       |
| 180                                  | DVI Cable      | SHIMADZU     | 2                | Yes      |
| 181                                  | DVI Cable      | SHIMADZU     | 2                | Yes      |
| 182                                  | DVI Cable      | SHIMADZU     | 2                | Yes      |
| 183                                  | DVI Cable      | SHIMADZU     | 2                | Yes      |
| 184                                  | AC Power Cable | SHIMADZU     | 2                | No       |
| 185                                  | AC Power Cable | SHIMADZU     | 2                | No       |
| 186                                  | AC Power Cable | SHIMADZU     | 2                | No       |
| 187                                  | AC Power Cable | SHIMADZU     | 2                | No       |

## ■ Guidance and Manufacturer's Declaration-Electromagnetic Emissions

| Guidance and Manufacturer's Declaration-Electromagnetic Emissions   |                |  |
|---|----------------|--|
| DAR-9500f is intended for use in the electromagnetic environment specified below. The customer or the user of DAR-9500f should assure that it is used in such an environment. |                |  |
| Emission Test   | Compliance     | Electromagnetic Environment-Guidance   |
| RF emissions<br>CISPR 11  | Group 1        | DAR-9500f uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                         |
| RF emissions<br>CISPR 11  | Class A        | DAR-9500f is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purpose. |
| Harmonic emissions<br>IEC 61000-3-2   | Not applicable |  |
| Voltage<br>fluctuations/flicker<br>emissions<br>IEC 61000-3-3   | Not applicable |  |


## Guidance and Manufacturer's Declaration-Electromagnetic Immunity

| Guidance and Manufacturer's Declaration-Electromagnetic Immunity  |   |  |  |
|---|---|--|--|
| DAR-9500f is intended for use in the electromagnetic environment specified below. The customer or the user of DAR-9500f should assure that it is used in such an environment. |   |  |  |
| Immunity Test   | IEC 60601 Test Level  | Compliance Level   | Electromagnetic Environment-Guidance   |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2  | ±6 kV contact<br>±8 kV air  | ±6 kV contact<br>±8 kV air   | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.   |
| Electrical fast transient/burst<br>IEC 61000-4-4  | ± 2kV for power supply lines<br>± 1kV for input/output lines  | ± 2kV for supply lines<br>± 1kV for input/output lines   | Mains power quality should be that of a typical commercial or hospital environment.  |
| Surge<br>IEC 61000-4-5  | ± 1kV line (S) to line (s)<br>± 2kV line (s) to earth   | ± 1kV line (S) to line (s)<br>± 2kV line (s) to earth  | Mains power quality should be that of a typical commercial or hospital environment.  |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11  | <5 % $U_T$<br>T(>95 % dip in $U_T$ ) for 0.5 cycle<br><br>40 % $U_T$<br>(60 % dip in $U_T$ ) for 5 cycle<br><br>70 % $U_T$<br>(30 % dip in $U_T$ ) for 25 cycle<br><br><5 % $U_T$<br>T(>95 % dip in $U_T$ ) for 5 s | Not applicable<br><br><br><br><br><br><br><br><br><br><5 % $U_T$<br>T(>95 % dip in $U_T$ ) for 5 s | Mains power quality should be that of a typical commercial or hospital environment. If the user of the DAR-9500f requires continued operation during power mains interruptions, it is recommended that DAR-9500f be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60Hz) magnetic field<br>IEC 61000-4-8   | 3 A/m   | 3 A/m  | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  |



Remark)  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and Manufacturer's Declaration-Electromagnetic Immunity

| <b>Guidance and Manufacturer's Declaration-Electromagnetic Immunity</b><br><b>DAR-9500f is intended for use in the electromagnetic environment specified below. The customer or the user of DAR-9500f should assure that it is used in such an environment.</b> |   |   |  |
|---|---|---|--|
| <b>Immunity Test</b>  | <b>IEC 60601 Test Level</b>                             | <b>Compliance Level</b>   | <b>Electromagnetic Environment-Guidance</b>  |
| Conducted RF<br>IEC 61000-4-6<br><br>Radiated RF<br>IEC 61000-4-3   | 3 Vrms<br>150 kHz 80 MHz<br><br>3 V/m<br>80 MHz 2.5 GHz | 3 Vrms<br>150 kHz to 1000 MHz<br><br>3 V/m<br>351.2 MHz<br>1980 MHz<br>2412 MHz | <p>Portable and mobile RF communications equipment should be used no closer to any part of DAR-9500f, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance:</b><br/> <math>d = 1.2\sqrt{P}</math></p> <p><math>d = 1.2\sqrt{P}</math> 80 MHz to 800 MHz<br/> <math>d = 2.3\sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |



At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.





These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which DAR-9500f is used exceeds the applicable RF compliance level above, DAR-9500f should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating DAR-9500f.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.



## WARNING

An exemption has been used any DAR-9500f has not been tested for radiated RF IMMUNITY over the entire frequency range 80 MHz to 2.5 GHz.

DAR-9500f has been tested for radiated RF IMMUNITY only at selected frequencies.

### List of the transmitters or equipment used as RF test sources and the frequency and modulation characteristics of each source.

| Kind of Equipment     | Type     | Manufacturer | Spot Check Frequency | Modulation  |
|-----------------------|----------|--------------|----------------------|---|
| For Radiated Immunity |          |              |                      |   |
| Digital Transceiver   | IC-DPR5  | ICOM         | 351.2MHz             | FSK (Frequency Shift Modulation)                  |
| Cellular Telephone    | 812SH    | Sharp        | 1980MHz              | PM (Phase Modulation)                             |
| Wireless LAN Station  | WHR-HP-G | BUFALLO      | 2412MHz              | OFDM (Orthogonal Frequency-Division Multiplexing) |

## Recommended Separation Distance Between Portable and Mobile RF Communication Equipment and the EQUIPMENT or SYSTEM

| Recommended Separation Distances Between Portable and Mobile RF Communication Equipment and DAR-9500f   |  |  |   |
|---|--|--|---|
| DAR-9500f is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of DAR-9500f can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and DAR-9500f as recommended below, according to the maximum output power of the communications equipment. |  |  |   |
| Rated Maximum Output Power of Transmitter<br>W  | Separation Distance According to Frequency of Transmitter<br>m |  |   |
|   | 150 kHz to 80 MHz<br>$d = 1.2\sqrt{P}$                         | 80 MHz to 800 MHz<br>$d = 1.2\sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 2.3\sqrt{P}$ |
| 0.01  | 0.12   | 0.12                                   | 0.23                                    |
| 0.1   | 0.38   | 0.38                                   | 0.73                                    |
| 1   | 1.2  | 1.2                                    | 2.3                                     |
| 10  | 3.8  | 3.8                                    | 7.3                                     |
| 100   | 12   | 12                                     | 23                                      |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.  |  |  |   |



At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.



These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**WARNING**

When using ACCESSORIES such as transducers and cables other than those specified, with the exception of those accessories sold by the Shimadzu of DAR-9500f as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of DAR-9500f.

**WARNING**

Do not use DAR-9500f next to or stack with other equipment. However, if it has to be used in that condition, DAR-9500f should be observed to verify that it operates normally.

### ■ Accessory List

The accessory list which conforms to EMC is shown below.

| Item                     | Model               | Manufacturer | Connection to |
|--------------------------|---------------------|--------------|---------------|
| Monitor                  | SMD 18101 SCS       | EIZO         | DAR-9500f     |
| Monitor                  | SMD 19102 SC        | EIZO         | DAR-9500f     |
| Monitor                  | RadiForce G11-S     | EIZO         | DAR-9500f     |
| Monitor                  | MVGD1519            | Barco        | DAR-9500f     |
| Monitor                  | ME181L/r-B          | TOTOKU       | DAR-9500f     |
| Monitor                  | RX150               | EIZO         | DAR-9500f     |
| Monitor                  | CL19194             | JVC KENWOOD  | DAR-9500f     |
| Monitor                  | ML19001             | JVC KENWOOD  | DAR-9500f     |
| Monitor                  | LS560W              | EIZO         | DAR-9500f     |
| Monitor                  | LX600W              | EIZO         | DAR-9500f     |
| Monitor                  | LS580W              | EIZO         | DAR-9500f     |
| Monitor                  | FlexScan EV2335W    | EIZO         | DAR-9500f     |
| Monitor                  | FlexScan S2100      | EIZO         | DAR-9500f     |
| Monitor                  | FlexScan S2133H     | EIZO         | DAR-9500f     |
| Monitor                  | S1934               | EIZO         | DAR-9500f     |
| Monitor                  | MX192               | EIZO         | DAR-9500f     |
| Monitor                  | CL19196             | JVC KENWOOD  | DAR-9500f     |
| Monitor                  | RX440               | EIZO         | DAR-9500f     |
| Monitor                  | RX850               | EIZO         | DAR-9500f     |
| Monitor                  | LCD-EA 193Mi-BM     | NEC          | DAR-9500f     |
| Monitor                  | MDSC-8258           | Barco        | DAR-9500f     |
| Large Monitor Manager    | LMM56800            | EIZO         | DAR-9500f     |
| Large Monitor Manager    | LMM0802             | EIZO         | DAR-9500f     |
| Large Monitor Controller | CID1000P            | EIZO         | DAR-9500f     |
| Large Monitor Controller | CID1201P            | EIZO         | DAR-9500f     |
| PC                       | Image Processing PC | HPC SYSTEMS  | DAR-9500f     |
| PC                       | GX620               | DELL         | DAR-9500f     |
| PC                       | Precision380        | DELL         | DAR-9500f     |
| PC                       | Precision390        | DELL         | DAR-9500f     |
| PC                       | Precision670        | DELL         | DAR-9500f     |
| PC                       | Precision690        | DELL         | DAR-9500f     |

| Item                    | Model            | Manufacturer       | Connection to |
|-------------------------|------------------|--------------------|---------------|
| PC                      | PrecisionT3400   | DELL               | DAR-9500f     |
| PC                      | PrecisionT5400   | DELL               | DAR-9500f     |
| PC                      | PrecisionT7400   | DELL               | DAR-9500f     |
| PC                      | PrecisionT3500   | DELL               | DAR-9500f     |
| PC                      | PrecisionT5500   | DELL               | DAR-9500f     |
| PC                      | PrecisionT7500   | DELL               | DAR-9500f     |
| PC                      | CELSIUS          | Fujitsu            | DAR-9500f     |
| PC                      | Express          | NEC                | DAR-9500f     |
| PC                      | PC, 3D-ANGIO     | HPC SYSTEMS        | DAR-9500f     |
| PC                      | GATEWAY2         | HPC SYSTEMS        | DAR-9500f     |
| PC                      | AR2200           | PFU                | DAR-9500f     |
| PC                      | AR8300           | PFU                | DAR-9500f     |
| PC                      | Express5800/56Xg | NEC                | DAR-9500f     |
| Local Image Controller  | IVR NEO          | SHIMADZU           | DAR-9500f     |
| Remote Image Controller | IVR Shuttle      | SHIMADZU           | DAR-9500f     |
| Dosimeter               | DIAMENTOR K2S    | PTW                | DAR-9500f     |
| Dosimeter               | DAP Meter        | SHIMADZU           | DAR-9500f     |
| Dosimeter               | VacuDAP          | VACUTEC            | DAR-9500f     |
| System Display          | System Display   | SHIMADZU           | DAR-9500f     |
| Intercom                | Interphone       | SHIMADZU           | DAR-9500f     |
| KVM Switch              | CS82A            | ATEN INTERNATIONAL | DAR-9500f     |
| KVM Switch              | CS1782A          | ATEN INTERNATIONAL | DAR-9500f     |
| SMART Touch             | SMART Touch      | SHIMADZU           | DAR-9500f     |

### ■ Essential Performance

- Image Acquisition
- Image Display
- Image Processing

### 2.6.7 Statement of Compliance (for Europe)

---

#### ■ Regulatory Information

For Europe:

The product complies with the requirement of the Medical Device Directive 93/42/EEC and RoHS Directive 2011/65/EU.

| Product Name                    | Digital Angiography System                                      |
|---------------------------------|---|
| Model Name                      | DAR-9500f   |
| Parts Number                    | 562-28000   |
| Manufacturer                    | SHIMADZU CORPORATION<br>Medical Systems Division                |
| Address                         | 1, NISHINOKYO-KUWABARACHO, NAKAGYO-KU, KYOTO<br>604-8511, JAPAN |
| Authorized Representative in EU | SHIMADZU EUROPA GmbH  |
| Address                         | Albert-Hahn-strasse 6-10, D-47269<br>Duisburg, F.R.Germany      |

#### ■ Company's Quality System

The company's Quality System is satisfied with Annex II, Article 3 for 93/42/EEC, which is certified by TUV Rheinland Product Safety GbmH (Notified under No.0197) as Registration No.: HD 60011592 001.

### 2.6.8 Other

---



The system needs activation of X-ray generator during the system installation.

If hospital name in the digital system is changed after the installation, activation is needed as well. (After changing hospital name, there is a certain grace period until activation.)

The activation work needs authorization by Shimadzu.

\*This is available for Trinias system shipped after February, 2019.

This system is connected to the hospital network via a router with firewall. The firewall function provides protection against threats to network security such as unauthorized access and viruses.

# Chapter3

## System Startup and Shutdown

This chapter describes how to start up and shut down the system.


### ■ Description of Chapter

|     |                          |      |
|-----|--------------------------|------|
| 3.1 | Startup .....            | 3-2  |
| 3.2 | User Authentication..... | 3-4  |
| 3.3 | Shutdown .....           | 3-10 |
| 3.4 | Power Failure .....      | 3-13 |

## 3.1 Startup


Follow this procedure to start up the system.

- 1 Press the system display ON button.



All the devices including the Digital processing system and the C-arm are activated.


- 2 As the GUI starts up, it briefly displays the software splash screen which includes information such as software version.



If the system ended while displaying a study, the study restarts automatically.

- 3 If using a Temporary License, the Temporary License status window is displayed. The window will be closed if fluoroscopy or new study starts.



If you do not want to use a temporary license and the display bothers you, enable to expire it immediately on the license administration window displayed when pressing [License Info] button.

Refer to  "15 Temporary License" P.15-1 for details.



**NOTE**

If installation of X-ray generator is not completed or the hospital name has changed, the following message may be displayed. If a message is displayed, please contact our service representatives.



## 3.2 User Authentication

If the Cyber Security feature is enabled, you can manage user accounts and record audit logs.



**NOTE**

This is initially disabled. Please contact our service representative to enable it.

### 3.2.1 User Level

There are two user levels: "USER" for users and "SUPER" for administrators. To switch user levels, sign in using the user ID and password for each user level.

If you sign in with an administrator user ID, you can configure cybersecurity settings.


### 3.2.2 Sign In

You must sign in to use the system.


Sign-in is required when you start a new study or using the Studies Management window.

[At the start of New Study]

- 1 Click [New] from the side menu.




- 2 The Sign In window appears.




[When using Studies Management window]


- 1 Click [Open] from the side menu.



- 2 Click the Sign In button on the Studies Management window.




- 3 The Sign In window appears.



### 3.2.3 Input User ID and Password


- 1 On the Sign In window, enter the desired user level "User ID" and "Password".



If you forgot your User ID or Password, click "Forgot User ID or Password?".  
You can sign in with user rights by entering the name and reason.

👉 "3.2.6 Emergency Sign In" P.3-8


- 2 Change the password only when you sign in for the first time. Enter "New Password" and "Confirm Password".



The image shows a 'Sign In' dialog box with a dark background. It contains two text input fields: 'New Password' and 'Confirm Password'. Below the fields are two buttons: 'OK' and 'Cancel'.



If the sign-in failure reaches the "Password failure times for user account locking" configured in the password policy, please restart the application.



The image shows a 'Sign In' dialog box with a dark background. It contains two text input fields: 'User ID' (containing 'guestuser') and 'Password'. Below the fields is a message: 'Exceeded failure time of password. Please restart the application.' At the bottom left is a link: 'Forgot User ID or Password?'. At the bottom right are two buttons: 'OK' and 'Cancel'.


[At the start of New Study]

After sign in, the New Study window appears.

[When using Studies Management window]

After sign in, the Study Management window is enabled.

In addition, [Sign In] button changes to [Sign Out] button and displays a new User ID.



The image shows the 'Studies Management' window. At the top, there is a header bar with the title 'Studies Management' and a user ID 'guestuser'. Below the header, there is a toolbar with buttons for 'System', 'Drives', 'Search', 'Change Pwd.', 'Sign Out', and 'Close'. The main area of the window contains a table with columns: 'Accession Number', 'Patient Name', 'Patient ID', 'Physician Name', 'Date', 'Origin', and 'Modality'. The 'Date' column has a value '02/19/2024' and a dropdown arrow. Below the table, there is a large empty space.

### 3.2.4 Change Password

Enable to change password.



**NOTE**

If the password change prohibition period is set, you cannot change the password during that period.

**Change Password**

Old Password


New Password

Confirm Password

Changing your password is currently prohibited. See the policies menu.

OK Cancel

- 1 Click [Change Pwd.] button.



- 2 Change the Password window appears. Enter "Old Password", "New Password", and "Confirm Password".

**Change Password**

Old Password

New Password

Confirm Password

OK Cancel

### 3.2.5 Sign Out

---

Sign out when you are finished with the system running.

Click [Sign Out] button on the Studies Management window.



If it is out of study, you will be signed out if you exit the system while signed in.


If the system restarts during study, hold the sign-in status.

### 3.2.6 Emergency Sign In

---

If you forgot your User ID or Password, you can sign in with user rights by entering the name and reason. Emergency Sign In is intended for temporary use when it is impossible to contact the administrator. Please make sure to have the password changed from the administrator account.


- 1 Click [Sing In].
- 2 Click [Forget User ID or Password?].

A screenshot of a 'Sign In' dialog box. It has fields for 'User ID' and 'Password'. At the bottom, there is a link 'Forgot User ID or Password?' which is highlighted with a red rectangular border. To the right of the link are 'OK' and 'Cancel' buttons.

- 3 Enter "Full Name" and "Reason", then clicc [OK].

A screenshot of an 'Emergency Sign In' dialog box. It has a 'Full Name' field and a larger 'Reason' text area. Both fields are highlighted with blue rectangular borders. At the bottom right are 'OK' and 'Cancel' buttons.

4 Then the user account name is displayed.



3




The user level for "emergency" is "USER". User account management features is not available.


# 3.3 Shutdown


Follow this procedure to shut down the system:

- 1 On the Reference monitor side menu, click [Exit].




- 2 If there are images in the print list that have not yet been printed, you can choose to print them now.




Click [No], to discard the images from the print list. Click [Yes] to print the images. Print window is displayed. Refer to  "Print Images" P.7-21.

- 3 If a network transfer or CD/DVD read is in progress, you can wait for it to finish or cancel it. Click [No] to allow it to finish. Click [Yes] to cancel it; loops and images already transferred, will be available at next startup.




- 4 If media (CD/DVD) write is in progress, a message appears indicating that you can either wait for completion or cancel it. Click [OK] to clear the message. If possible, wait for the media write to finish.






- 5 If necessary, you can click [Cancel Media Write] on the side menu.




When prompted, click [Yes] to cancel the write and discard the media or click [No] and wait for the write to finish.




If the installation personnel or administrator configured a shell command to be executed upon shutdown, and the command fails, a message to this effect will remain on screen until it is acknowledged.

- 6 If a study is open for acquisition, a prompt appears. Click [Yes] to close the study now. Otherwise, click [No] to leave the study open (study will be continued at next startup).



- 7 Finally, you are prompted to complete the shutdown. Click [Yes] to complete the shutdown.



- 8 Wait for the Digital processing system to be completely off (display monitors go blank and green power indicators on cabinets go off).




#### NOTE


If trouble occurs when trying to shut down the Digital processing system via the GUI [Exit] button (system does not shut down), press and hold the Digital reference cabinet power button for at least five seconds until the power goes off (green indicator goes off). Press [ACQ POWER F] button of Digital acquisition cabinet in the Equipment room to turn off the power (indicated as ACQ PC in the cabinet).

For Bi-plane system, also press [ACQ POWER L] button of Digital acquisition cabinet to turn off the power.

To restart the system, wait at least one minute after turning the power off the System Display, and press the ON button again.



## 9 Press the System display [OFF] button.



#### NOTE

If something occurs to the application and reboots the DAR-9500f application, press the [Tab] key while holding the [Alt] key on the keyboard. Then select [System Reboot], and click [OK] two times. The DAR-9500f application reboots.

## 3.4 Power Failure

When a power failure occurs, the GUI disappears. If the power is restored within 5 minutes, the GUI is displayed again and is immediately functional. However, if the power failure lasts more than 5 minutes, wait at least 10 minutes and startup the system as usual.



This system is designed to minimize any data loss due to power failure. Only the specific image loop that was in acquisition at the time of the power failure is lost. If the power is restored after more than 10 minutes of power failure, the stud which was opened for acquisition right before the power failure restarts.



If an instant power failure or a power failure occur, FPD calibration returns to default value and an appropriate image may not acquire until auto calibration is started.

Perform calibration manually after recovering from the power failure.

"19.3.1 Calibration" P.19-5



Shimadzu recommends to contact a service personnel when a power failure occurs.



Startup of the system within 10 minutes of a power restore may damage the system.

This page is intentionally left blank.

# Chapter 4

## Image Acquisition

This chapter introduces the main user interface window and describes how to acquire loops and images during a medical diagnostic procedure as follows:

### Description of Chapter

|      |  |       |
|------|--|-------|
| 4.1  | Introduction .....                                 | 4-2   |
| 4.2  | The Main User Interface Window. ....               | 4-2   |
| 4.3  | Entering Study Information .....                   | 4-10  |
| 4.4  | Making Acquisition. ....                           | 4-23  |
| 4.5  | Actions When Study Open for Acquisition .....      | 4-50  |
| 4.6  | Working With Reference Images .....                | 4-60  |
| 4.7  | Acquisition Keys and GUI Buttons .....             | 4-63  |
| 4.8  | IVR NEO (Either-or SMART Touch) .....              | 4-65  |
| 4.9  | Keyboard ShortCuts. ....                           | 4-72  |
| 4.10 | IVR Shuttle (Option) (Either-or SMART Touch) ..... | 4-73  |
| 4.11 | SMART Touch (Either-or IVR NEO/IVR Shuttle) .....  | 4-77  |
| 4.12 | Closing the Active Study .....                     | 4-105 |


## 4.1 Introduction

DAR-9500f acquires loops during X-ray angiography diagnostic procedure. As a loop is acquired, it appears on the Acquisition monitor and soon after on the Reference monitor. Rad loops are recorded and Fluoro loops may optionally be recorded. A use-selected reference image can be selected and included with each loop. Loops from previous studies for the current patient can be displayed on the Reference monitor for comparison against the current loops displayed on the Acquisition monitor.


## 4.2 The Main User Interface Window

The main DAR-9500f graphical user interface (GUI) window is displayed on the Reference monitor and the Acquisition monitor. It features a Side Menu, an Image Viewer, an Image Selector, and a Cine Control and Status bar along the bottom. Additional elements such as Digital Program Pop-up Menus and acquisition parameters are superimposed over the Image Viewer window.

- Acquisition Monitor Main Window



- Reference Monitor Main Window







The above sample screens are made from a composite of images. Not all elements appear at the same time.


## 4 Image Acquisition




The key elements of the main DAR-9500f user interface window, numbered individually for the Reference and Acquisition monitors in the above illustrations, are defined in the following table.












**A:** Display on Acquisition monitor, **R:** Display on Reference monitor.

| No. | Item                | Description  | A                     | R                     |
|-----|---------------------|--|-----------------------|-----------------------|
| 1   | Side Menu           | Provides access to all major functionality.  | <input type="radio"/> | <input type="radio"/> |
|     | Reference Indicator | (Not shown) When a monitor is in Reference mode (IVR NOE/IVR Shuttle), the indicator "Ref" appears in its upper-left corner.   |                       |                       |
| 2   | Study Information   | <p>Displays study information as follows:</p> <ul style="list-style-type: none"> <li>Image/Loop, expressed as: FRAME/TOTAL_FR-IMAGE[TOTAL_IM] (IM_TYPE)<br/> FRAME= Frame number<br/> TOTAL_FR= Total number of frames in loop<br/> IMAGE= Current image/loop number<br/> (A for plane A or B for plane B)<br/> TOTAL_IM= Total number of loops<br/> IM_TYPE= Type of images<br/> (Such as DA, DSA, 3D-DA)<br/> (LL= Lossless (if not LL, then raw)</li> <li>Study Accession Number</li> <li>Patient Name</li> <li>Patient Age</li> <li>Patient Sex</li> <li>Patient Birth Date</li> <li>Patient ID</li> <li>Study ID</li> <li>Study Date and Time</li> </ul> <p> <a href="#">"4.3 Entering Study Information" P.4-10</a></p> | <input type="radio"/> | <input type="radio"/> |
| 3   | Acquisition Active  | <p>Indicates the active acquisition type such as: Fluo, MAP, DA, DSA, RSM</p> <p> <a href="#">"4.4.4 Making Live Fluoro Acquisitions" P.4-30</a></p> <p> <a href="#">"4.4.6 Recording Rad" P.4-40</a></p>  | <input type="radio"/> |                       |
| 4   | Image Selector      | <p>Displays one icon per loop or still image. For loops, the icon is automatically created from the image at one-third of the way into the loop. For example, if loop consisting of 270 images, the icon is created from the 90 th image.</p> <p> <a href="#">"5.3 Displaying Study Information" P.5-6</a></p>  | <input type="radio"/> | <input type="radio"/> |



| No. | Item       | Description   | A                     | R                     |
|-----|------------|---|-----------------------|-----------------------|
| 5   | Parameters | Dynamically-updated Information (individually configured):  |                       |                       |
|     |            | mGym <sup>2</sup> :<br>Cumulative X-ray dosage since the beginning of study (Acquisition monitor).<br>Displayed loop X-ray dosage (Reference monitor).  | <input type="radio"/> | <input type="radio"/> |
|     |            | mGy/min:<br>(AKR)(Air kerma rate) Dose per unit of time.<br>The total X-ray dosage (Rad and Fluoro) at Patient Entrance Reference Point (Interventional Reference Point), that would occur if the active exposure was continued for a minute. AKR is displayed only while performing fluoroscopy and radiography.                         | <input type="radio"/> |                       |
|     |            | mGy:<br>(AK)(cumulative air kerma) The cumulative dose of all X-ray exposures (Rad and Fluoro), since start of study or manual dose reset at Patient Entrance Reference Point (Interventional Reference Point).<br>"Show dose warning" has been set, and if the dose exceeds the threshold, the underline and warning icon are displayed. | <input type="radio"/> |                       |
|     |            |    |                       |                       |
|     |            | kV:<br>X-ray tube voltage in kilo-volts.  | <input type="radio"/> | <input type="radio"/> |
|     |            | mA:<br>X-ray tube current in milliamperes.  | <input type="radio"/> | <input type="radio"/> |
|     |            | ms:<br>X-ray pulse duration per image, in milliseconds.   | <input type="radio"/> | <input type="radio"/> |
|     |            | I-A (s):<br>The time in seconds between injection of contrast agent and start of acquisition.   | <input type="radio"/> | <input type="radio"/> |
|     |            | AcqTime:<br>Acquisition time.   | <input type="radio"/> | <input type="radio"/> |
|     |            | LAO-RAO:<br>Oblique angle of the detector in degrees in which angles toward the patient's left are displayed with LAO, and angles toward the patient's right are displayed with RAO.  | <input type="radio"/> | <input type="radio"/> |
|     |            | CRA-CAU:<br>Sagittal angle of the detector in degrees in which angles toward the patient's head are displayed with CRA, and angles toward the patient's foot are displayed with CAU.  | <input type="radio"/> | <input type="radio"/> |
|     |            | SID (cm):<br>X-ray source to detector distance in centimeters.  | <input type="radio"/> | <input type="radio"/> |

| No. | Item                 | Description   | A                     | R                     |
|-----|----------------------|---|-----------------------|-----------------------|
|     |                      | <p>FOV (in):<br/>Detector field of view in inches (diameter for image intensifiers).</p> <p>HDD (min):<br/>Possible remaining radiography time based on the hard disk capacity.</p> <p>Fluo:<br/>The total number of minutes of Fluoro exposure since start of study or manual reset via generator console.</p> <p>I-F:<br/>The period of time when turn ON the injector with overlay by frame.</p> <p>I-T (s):<br/>Cumulative time in seconds since the beginning of injection of contrast agent.</p> <p>&lt;seconds&gt;:<br/>(not shown) (at bottom right) FluoroMAP Timer (seconds) counts down at the beginning of a FluoroMAP acquisition.</p> <p>LIVE/LIH/STORED:<br/>The status of an image.</p> <p>Patient Direction:<br/>Direction of a patient body is indicated as follows:<br/>H=Head<br/>F=Feet<br/>R=Right<br/>L=Left<br/>A=Anterior<br/>P=Posterior<br/> <a href="#">"Display Tab" P.17-37</a></p> <div style="border: 1px solid black; padding: 10px; margin-top: 10px;"> <div style="text-align: center;">  <b>NOTE</b> </div> <p>When acquiring the image by rotating the C-arm, the C-arm angle can be saved in each frame. However, there is an error up to 3 degrees more or less between angles to be saved in frame and actual C-arm angle depending on the system.</p> <p>If save or display angles after moving the C-arm, please contact our service representative.</p> </div> | <input type="radio"/> | <input type="radio"/> |
| 6   | Digital Program Menu | <p>Selects Digital Programs for acquisition.</p> <p> <a href="#">"4.4.3 Choosing Acquisition Modes" P.4-26</a></p>   | <input type="radio"/> | <input type="radio"/> |
| 7   | Probe Signal         | <p>During review, up to four probe signals, including ECG and blood pressure, can be displayed as curves superimposed over the lower half of the Image Viewer window. A moving cursor (vertical line) is synchronized to the currently-displayed frame.</p>   | <input type="radio"/> | <input type="radio"/> |

| No. | Item            | Description   | A                     | R                     |
|-----|-----------------|---|-----------------------|-----------------------|
| 8   | Image Viewer    | The main area for viewing and working with images.<br> <a href="#">"4.2 The Main User Interface Window" P.4-2</a>  | <input type="radio"/> | <input type="radio"/> |
| 9   | Timer           | The stopwatch-like timer displays from 0 to 59 min with a resolution of 0.1 s. Click the [Timer] side-menu button once or the IVR NEO/IVR Shuttle/SMART Touch [Timer] button once to start the timer, again to stop it, and a third time to hide it.  | <input type="radio"/> | <input type="radio"/> |
| 10  | Status Messages | Displayed when needed in both the Image Viewer and Studies Management windows.<br> <a href="#">"6.1.7 Messages &amp; Information" P.6-8</a>  | <input type="radio"/> | <input type="radio"/> |
| 11  | Programs        | Displays names of selected Rad and Fluoro digital programs.<br> <a href="#">"4.4.3 Choosing Acquisition Modes" P.4-26</a>  | <input type="radio"/> | <input type="radio"/> |
| 12  | Exit Button     | Exit the system.  | <input type="radio"/> | <input type="radio"/> |
| 13  | Cine Control    | Provides transport controls for loop playback.<br><br> <a href="#">"4.7 Acquisition Keys and GUI Buttons" P.4-63</a><br> <a href="#">"5.2.1 Controlling Image Playback" P.5-5</a><br> <a href="#">"7.1.1 Mouse and Keyboard" P.7-2</a>  | <input type="radio"/> | <input type="radio"/> |
| 14  | Modes & Info.   | These buttons appear as needed in this order: Display Study Information, Set Playback Mode, and Display Error List;<br><br>Click  to display the list like this.<br><br>Available buttons are shown as follows.<br>[Manage Error]: Enable when selecting an error of MPPS support.<br>Click [Manage Error] to cancel retry of MPPS support.<br>[Copy]: Select any message to make a text copy of the message.<br>[Delete]: Select any message to delete it.<br>[Delete All]: Delete every message.<br>[Close]: Close the error list.<br> <a href="#">"6.1.7 Messages &amp; Information" P.6-8</a> | <input type="radio"/> | <input type="radio"/> |




| No. | Item             | Description   | A                     | R                     |
|-----|------------------|---|-----------------------|-----------------------|
| 15  | Network Activity | <p>Displays a simple animation during network activity. Also, click this area to display a Network Status box.</p>  <p>To hide the Network Status box, click the Network Activity area (lower right of bottom bar) again.</p> <p> <a href="#">"6.1.7 Messages &amp; Information" P.6-8</a></p> | <input type="radio"/> | <input type="radio"/> |
| 16  | Comment          | <p>The Comment (if any) selected from the list of site-specific comments. (When you select the comment in advance, you can acquire a loop with this comment displayed for the next acquisition.)</p> <p> <a href="#">"4.4.2 Adding a Comment" P.4-25</a></p>   | <input type="radio"/> | <input type="radio"/> |



Image edges may be blacked out by shutters.

### 4.2.1 Reference and Acquisition Monitors

Most GUI interaction occurs on the Reference monitor.

The keyboard and mouse normally control the Reference monitor.

 ["4.7 Acquisition Keys and GUI Buttons" P.4-63](#)

Some keyboard actions can be performed on the Acquisition monitor by holding down the Keyboard [Ctrl] key at the same time as pressing another keyboard key.

 ["7.1.1 Mouse and Keyboard" P.7-2](#)

The IVR NEO/IVR Shuttle/SMART Touch [Monitor Select] button enables the Reference or Acquisition monitor to be chosen. IVR NEO/SMART Touch joystick and button actions apply to the selected monitor.

 ["4.8 IVR NEO \(Either-or SMART Touch\)" P.4-65](#), ["4.10 IVR Shuttle \(Option\) \(Either-or SMART Touch\)" P.4-73](#), ["4.11 SMART Touch \(Either-or IVR NEO/IVR Shuttle\)" P.4-77](#)



Unless otherwise stated, all GUI-interaction descriptions in this manual are for the Reference monitor.

## 4.3 Entering Study Information

### 4.3.1 Introduction


---

A new study must be defined before image acquisition can be performed.


### 4.3.2 Defining a New Study

---

To define a new study, click [New] on the side menu.




If installation of X-ray generator is not completed or the hospital name has changed, the following message may be displayed. If a message is displayed, please contact our service representatives.



**Define New Study with Worklist Server Management**

If the [Worklist Server Management] window appears, your institution uses a worklist server system that automates much of the new-study-definition process.

If an empty [New Study] window appears, continue with  ["Define New Study with Direct Entry"](#) [P.4-17](#).

WorkList Management

Search Scheduled Procedure Steps

Scheduled Procedure Step Start Date:

Month

Day

Year

Month

Day

Year

From

02

14

2018

To

02

14

2018

Modality Worklist Server

MWMSCP

Modality

XA

Patient Name

Requested Procedure ID

Patient ID

Station AE Title

Performing Physician's Name

Accession Number

Find

Cancel Find

☐ Auto Refresh

☐ Auto Fill AE Title

List of Scheduled Procedure Steps

Number of Scheduled Procedure Steps found: 0

| Accession Number | Patient Name | Patient ID |
|------------------|--------------|------------|
|------------------|--------------|------------|

Select From MWL

Direct entry

Cancel

Follow this procedure to define a new study using worklist server management.

- 1 In the [Worklist Management] window, select the [Modality Worklist server] (if the default is not desired), enter an appropriate date range and any other search criteria and then click [Find].

After several seconds, [Procedure Steps] that match our criteria are displayed in the bottom half of the [Worklist Management] window. If check [Auto Fill AE Title], [DICOM AE Title] is set automatically to [Station AE Title]. If desired, click [Auto Refresh] (check mark appears) to cause automatic periodic refresh of the list.

The screenshot shows the 'WorkList Management' window. It has a search section at the top with fields for 'Modality Worklist Server' (set to 'MWMSCP'), 'Modality' (set to 'VA'), and 'Scheduled Procedure Step Start Date' (From: 02/26/2018, To: 02/26/2018). There are also fields for 'Patient Name', 'Requested Procedure ID', 'Patient ID', 'Station AE Title', 'Performing Physician's Name', and 'Accession Number'. Below these are 'Find' and 'Cancel Find' buttons, and checkboxes for 'Auto Refresh' and 'Auto Fill AE Title'. The bottom section, 'List of Scheduled Procedure Steps', shows 'Number of Scheduled Procedure Steps found: 4' and a table with 4 rows of data.

| Accession Number | Patient Name     | Patient ID |
|------------------|------------------|------------|
| 37               | SHIMADZU TARO037 | SHIMA0038  |
| 43               | SHIMADZU TARO043 | SHIMA0044  |
| 44               | SHIMADZU TARO044 | SHIMA0045  |
| 45               | SHIMADZU TARO045 | SHIMA0046  |

At the bottom of the window are buttons for 'Select From MWL', 'Details', 'Direct entry', and 'Cancel'.



### NOTE

If you cannot find the desired worklist server, click [Direct entry] instead of [Select From MWM], or just double-check the desired procedure step. Continue with ["Define New Study with Direct Entry" P.4-17](#).



### NOTE

If you click [Search] during media (CD/DVD) write, a message to stop media writing (discard the media) will be displayed before start searching. Stop writing the media, and click [Yes] to search and [No] to complete the media write.



- 2
- Look through the procedure steps and select one.  
It becomes highlighted in the list.
- 3
- If MPPS support is enabled, click [Details] to display any available information about procedure.

Scheduled Procedure Step Details

Modality:

KA

Station AE Title:

MWMSCU

Performing Physicians:

Performing PHYS

Scheduled Protocol Codes:

Scheduled Procedure Step Description:

Study Description

Requested Procedure Description:

RQP\_DESC

Requested Procedure Codes:

Requested Procedure ID:

RQP\_ID1

Referring Physician's Name:

REF\_PHYS4

Scheduled Procedure Step Start Date:

2018/02/26 16:16:26

Patient's Birth Date:

1965/12/26

Patient's Gender:

Female

Medical Alerts:

Close

- 4
- Click [Close].  
Make sure that the correct procedure step is selected in the list, and then click [Select From MWM], or just double-check the desired procedure step.
- 5
- The [New Study] window is displayed with all available information pre-entered.  
Items in gray cannot be changed. The [Free Space] of your hard disk is indicated at the top of the [New Study] window. (Without MPPS Support (left) and with (right))

New Study      Free Space: 34143MB

Patient

Patient ID

PIDDAR9500REV\_022

Generate

Month

07

Day

07

Year

1940

Age

79

Patient Name

Prefix

First Name

John

Middle Name

Last Name

Shimadzu

Suffix

Height

inches

Weight

lbs

Sex

Male

Select Orientation

Clear

Study

07/12/2019 - 14:14:38

Study ID

S00001

Accession Number

A00001

Modality

KA

Referring Physician

Referring First

Study Description

Performing Physician's List

☐ Daily Check  
☐ Performing First  
☐ Performing Second

Operators' List

☐ Operator First  
☒ Operator Second

Primary Physician: Default

☒ Auto Transfer    ☐ H Flip  
☒ Acquire Waves    ☐ V Flip

OK

Cancel

New Study      Free Space: 34135MB

Patient

Patient ID

PIDDAR9500REV\_022

Generate

Month

07

Day

07

Year

1940

Age

79

Patient Name

Prefix

First Name

John

Middle Name

Last Name

Shimadzu

Suffix

Height

inches

Weight

lbs

Sex

Male

Select Orientation

Clear

Study

20180831 105959.000

Study ID

10151081

Accession Number

55994001

Modality

KA

Referring Physician

Study Description

ScheduledProcedureStepDescription

Performing Physician's List

☐ Daily Check  
☐ Performing First  
☐ Performing Second

Operators' List

☐ Operator First  
☐ Operator Second


Primary Physician: Default


☒ Auto Transfer    ☐ H Flip  
☒ Acquire Waves    ☐ V Flip

OK

Cancel

- 6 If MPPS support is enabled, you can change the list of codes appearing in the Scheduled Protocol Codes list at the bottom of the New Study window.

Click the right-arrow at the right edge of the list, and then check/uncheck the codes as desired. Checked codes appear in the list; unchecked codes do not. The list of available codes is defined by the institution.  ["DICOM Tab" P.17-24](#)




Click [Modify] to save any changes, otherwise [Cancel]. The Scheduled Protocol Codes window closes and the New Study window re-appears.

- 7 If desired, choose options in available lists such as [Performing Physician], fill in any blank field such as [Protocol Name], and select options like [Auto Transfer] and [Acquire Waves], as appropriate. (Without MPPS support (left) and with (right))

The field marked with a dot in the upper corner are considered mandatory, and must be filled in by either the worklist server manager or you.

Multiple performing physicians can be selected. If so, a primary physician is selected as default. Right-click on the performing physician's name to change the primary physician.



If MPPS support is enabled, take Protocol Name off from mandatory items. In this case, DUP group name will be used for Protocol Name as default.

👉 "DICOM Tab" P.17-24

To automatically transfer studies to the default server, select [Auto Transfer]. Otherwise, ensure that [Auto Transfer] is deselected.

👉 "Physicians Tab" P.17-44

To acquire waves (ECG) during a case, select [Acquire Waves]. Otherwise, ensure that [Acquire Waves] is deselected.

👉 "Display Tab" P.17-37

- 8 Click [OK] when finished entering new study information.

If any mandatory fields have not been filled in by the worklist server manager or you, the [New Study] window opens in direct-entry mode enabling you to complete the mandatory fields. Complete such fields and click [OK]. The [New Study] window closes and the new study is opened and prepare for acquisition.

- 9 Continue with either  "Display New Study Information" P.4-21 or  "4.4 Making Acquisition" P.4-23.

## Define New Study with Local Worklist Management

Valid the local worklist to display the [Local Worklist Management] window.

- **Registration of New Study**

Registration of new study is required in advance when using the local worklist. Follow the instruction below to input new study data.

[illegible]


- 1 Select [Add Patient].

- 2 [Add New Patient] window is displayed. Input required items and select [Save].

- **Selection of New Study**

- 1 Click [New] on the side menu.
- 2 Patient information registered in the local study list is displayed. Registered patient can be searched by accession No. Patient Name, Patient ID and Date of Birth. Enter the keyword in the search at the top of window and click [Find].
- 3 Select the patient and click [Select From List].
- 4 [New Study] window is displayed with registered patient data. Gray-out items are automatically entered and cannot fix them. "Free Space" at the top of [New Study] window indicates the capacity of hard disk.

- 5
- Select [Performing Physician's List] and fill in the empty items if necessary.  
Multiple performing physicians can be selected. If so, a primary physician is selected as default. Right-click on the performing physician's name to change the primary physician.



- 6
- Click [OK] after input study information.
- 7
- Continue with either ["Display New Study Information" P.4-21](#) or ["Making Acquisition" P.4-23](#).

**Define New Study with Direct Entry**


When there is no worklist manager or you clicked the [Direct entry] button on the [Worklist Management] window, a blank [New Study] window appears. The "Free Space" of your hard disk is indicated at the top of the [New Study] window. (Without MPPS support (left) and with (right))




**NOTE**

When MPPS support is enabled and you choose [Direct Entry], the Accession Number field is not available (box background is gray).

With either blank [New Study] window, follow this procedure to define a new study without worklist management


- 1 If you worked with the patient before, you can select the patient from the Patient list. Click the right-arrow  to the right of the Suffix box. A patient list appears.




The dialog box titled "Patient" contains input fields for Last Name, First Name, and Patient ID, along with "Select" and "Close" buttons. Below these is a table of patient records.

| Last Name | First Name | Patient ID       | Birth Year | Birth Month | Birth Day | Sex |
|-----------|------------|------------------|------------|-------------|-----------|-----|
| Shimadzu  | Amy        | PIDDAR9500REV... | 1992       | 12          | 8         | F   |

- 2 If desired, reduce the number of patients shown by entering search criteria in the provided boxes: Last Name, First Name, or Patient ID. Only patients matching what you enter will be displayed.
- 3 If you see the desired patient, click anywhere in its row and then click [Select]. This loads all available patient information into the [New Study] window, reducing the number of fields you must enter manually. (Without MPPS support (left) and with (right))



The "New Study" window shows fields for Patient ID (PIDDAR9500REV\_022), Patient Name (John Shimadzu), Date of Birth (07/07/1940), Height, Weight, Sex (Male), and Study ID (500001). It includes checkboxes for "Performing First" and "Performing Second" under the "Performing Physician's List".



The "New Study" window shows fields for Patient ID (PIDDAR9500REV\_022), Patient Name (John Shimadzu), Date of Birth (07/07/1940), Height, Weight, Sex (Male), and Study ID (10151081). It includes checkboxes for "Performing First" and "Performing Second" under the "Performing Physician's List".


- 4 If the patient is not available in the list, you must manually enter the patient information in the [New Study] window.

Thereafter, the patient will be available in the patient list.

 "Patient Orientation" P.4-19

- 5 Choose [Referring Physician], [Performing Physician] and [Operator] from the available lists.

Multiple performing physicians can be selected. If so, a primary physician is selected as default. Right-click on the performing physician's name to change the primary physician.



The dialog box shows a list of performing physicians: "Performing First" and "Performing Second". A "Set Primary" button is next to "Performing First".

- 4
- 6

Fill in at least all blank mandatory fields (marked with a dot in upper corner).

For information on how the administrator can configure which fields are mandatory, see ["Study Information Tab" P.17-50](#).

If desired, click [Generate] to generate a unique Patient ID. If desired, click [Clear] to clear all fields in the Patient area of the dialog box.

To automatically transfer studies to the default server, select [Auto Transfer]. Otherwise, ensure that [Auto Transfer] is deselected.

["Physicians Tab" P.17-44](#)

To acquire waves (ECG) during case, select [Acquire Waves]. Otherwise, ensure that [Acquire Waves] is deselected.

["Display Tab" P.17-37](#)
- 7

When done, click [OK].

The [New Study] window closes and the new study is opened and prepared for acquisition.



If you choose a patient and then make changes to their information in the [New Study] window, the changes are retained in the patient record and will be seen in any future uses of the patient for new-study definition.

- 8

Continue with either the section ["Display New Study Information" P.4-21](#) or ["4.4 Making Acquisition" P.4-23](#).


**Patient Orientation**

Set patient orientation during study if necessary.


- 1

Click [Select Orientation] button.

- 2 Patient positions are listed as follows. Select an appropriate position.



- 3 Selected position is displayed as follows on [New Study] window. And H/V Flip is checked automatically depending on the selected position.


A screenshot of the 'New Study' window. The 'Patient' section shows fields for Patient ID (PIDDA9500REV\_022), Date of Birth (7/7/1940), Patient Name (John Shimadzu), Height, Weight, Sex (Male), and a position dropdown menu set to 'FFS'. The 'Study' section includes fields for Study ID, Accession Number, Modality (XA), Referring Physician, and Study Description. Below these are lists for 'Performing Physician's List' and 'Operators' List'. At the bottom, there are checkboxes for 'Auto Transfer', 'Acquire Waves', 'H Flip', and 'V Flip'. The 'H Flip' and 'V Flip' checkboxes are circled in red, indicating they are checked. 'OK' and 'Cancel' buttons are at the bottom right.

[Patient Orientation] can be selected only one position for one study. Do not select the position if changing patient orientation during study.



■ Display New Study Information

To display new study information, click [Study] on the side menu.



If MPPS support is enabled or any loops have been recorded or images have been saved or performed fluoroscopy for this study, the [Study Information] window appears. View the information and click [Close].

The screenshot shows a 'Study Information' window. It has sections for 'Patient' (with fields for ID, Name, Date of Birth, Height, Weight, Sex, Orientation), 'Study' (with fields for ID, Accession Number, Modality, Referring Physician), and 'Study Description'. There are also lists for 'Performing Physician's List' and 'Operators' List'. At the bottom, there are buttons for 'Anonymize', 'Modify', and 'Close'. The window title is 'Study Information'.

4

If MPPS support is not enabled or loops have not yet been recorded or images have not yet been saved/performed fluoroscopy, the [Study Modification] window is displayed instead. View the information, optionally make changed, and click [OK] to save changes or click [Cancel] to exit without saving. Items not available for modification are displayed in gray.



Completed studies can have much of their information modified as described in ["7.5 Modifying Study Information" P.7-12.](#)

### Study Modification

Patient

Patient ID

PIDDAR9500REV\_022

Generate

Date of Birth

☐ Unknown

Month

07

Day

07

Year

1940

Age

79

Patient Name

Prefix

First Name

Middle Name

Last Name

Suffix

John

Shimadzu

Select Orientation

Clear

Height

inches

Weight

lbs

Sex

Male

Study

07/12/2019 14:18:28

Study ID

S00001

Accession Number

A00001

Modality

KA

Referring Physician

Referring First

Study Description

Performing Physician's List

☐ Daily Check  
☒ Performing First  
☐ Performing Second

Operators' List

☐ Operator First  
☒ Operator Second

Primary Physician: Performing First

☒ Auto Transfer  
☒ Acquire Waves

☐ H Flip  
☐ V Flip

OK

Cancel

## 4.4 Making Acquisition




Any direction in this manual to expose Fluoro or Rad radiation are only for the purpose or procedure description. The operator must decide when it is medically safe to expose a patient to radiation.

### 4.4.1 Introduction


Fluoro acquisitions can be made and displayed live at any time (without recording). When a study is open for acquisition, all acquisition modes are available: **Fluoro, FluoroMap, SIMAP, DA, DSA, RSM-DSA, ROT-DA, ROT-DSA, ROT-RSM-DSA, SCORE 3D-DA, SCORE 3D-DSA, SCORE CT, HQ-DSA, Stage Acquisition and SPOT.**

If you have not already done so, define a new study as described in ["4.3.2 Defining a New Study" P.4-10](#). When a new study is defined, the [Image Viewer] window is blank except for a study information line at the top.

Display items of study information can be customized. See ["Study Information Tab" P.17-50](#).




Pediatrics setting is supported. The warning message will be displayed if patient is under the setting age. See ["Notifications Tab" P.17-42](#).













"Show dose warning" has been set, and if the dose exceeds the threshold, the underline and warning icon are displayed. See ["Notifications Tab" P.17-42](#).



For information about using the IVR NEO, IVR Shuttle and the GUI/Keyboard to assist with making acquisitions, see ["4.7 Acquisition Keys and GUI Buttons" P.4-63](#), ["4.8 IVR NEO \(Either-or SMART Touch\)" P.4-65](#), ["4.10 IVR Shuttle \(Option\) \(Either-or SMART Touch\)" P.4-73](#) and ["4.11 SMART Touch \(Either-or IVR NEO/IVR Shuttle\)" P.4-77](#). You can adjust image appearance of the recently-acquired loops as described in ["7.2 Adjusting Image Appearance" P.7-4](#).

The acquisition procedures are describe in the following sections.

| No. | Section                                 | Page                   |
|-----|---|------------------------|
| 1   | "4.4.2 Adding a Comment"                | <a href="#">P.4-26</a> |
| 2   | "4.4.3 Choosing Acquisition Modes"      | <a href="#">P.4-26</a> |
| 3   | "4.4.4 Making Live Fluoro Acquisitions" | <a href="#">P.4-30</a> |
|     | "FluoroMAP"                             | <a href="#">P.4-31</a> |
|     | "DSA-MAP"                               | <a href="#">P.4-32</a> |
|     | "BlankMAP"                              | <a href="#">P.4-34</a> |
|     | "LIVE Mode"                             | <a href="#">P.4-37</a> |
|     | "TraceMAP"                              | <a href="#">P.4-37</a> |
| 4   | "4.4.5 Recording Fluoro"                | <a href="#">P.4-38</a> |
|     | "Direct Fluoro Record"                  | <a href="#">P.4-38</a> |
|     | "Last N Seconds Fluoro Record"          | <a href="#">P.4-39</a> |

| No. | Section                  | Page   |
|-----|--------------------------|--|
| 5   | "4.4.6 Recording Rad"    |  P.4-40 |
|     | "DA"                     |  P.4-40 |
|     | "DSA"                    |  P.4-41 |
|     | "RSM-DSA"                |  P.4-42 |
|     | "ROT-DA and ROT-RSM-DSA" |  P.4-42 |
|     | "ROT-DSA"                |  P.4-43 |
|     | "SPOT"                   |  P.4-45 |
|     | "Flipping Images"        |  P.4-47 |

#### 4.4.2 Adding a Comment



To add comments, the comment menu must be configured.

See  ["Configure the Comment Menu" P.17-70](#).

You can add a predefined comment to radiography and recorded fluoroscopy acquisitions.

To add a comment, follow this procedure.

- 1 Click [Study]-[Comm.] GUI button or press the equivalent IVR NEO or SMART Touch button.



The comment menu appears, in small font if launched from the GUI buttons or in large font if launched from the IVR NEO.



- 2 Once a comment menu launches, use the keyboard or IVR NEO/SMART Touch joystick menu selector to make comment menu choices as follows.

| To move the selection highlight |                                     |
|---------------------------------|-------------------------------------|
| Keyboard:                       | Use the Up/Down arrow keys.         |
| IVR NEO/SMART Touch:            | Use the Up/Down motion of joystick. |

| To select a menu item |                            |
|-----------------------|----------------------------|
| Keyboard:             | Press [Enter].             |
| IVR NEO/SMART Touch:  | Press the joystick button. |

The highlighted menu item appears in black text with a white outline.

After selection, the menu is closed. Your selection in the menu is remembered until the study is closed, so that the same comment remains selected for all subsequent loops. When you select a comment, it is displayed at the bottom of the [Image Viewer] window as shown, for radiography and fluoroscopy acquisitions, playback and review, according to your specific settings.


See  "Display Tab" P.17-37.


Cardio / CAG-F[15f-10s]  
Cardio / 15pps  
Comment no.1

- 3 Alternatively, select [Cancel] or press the [End] key (keyboard) to close the menu without choosing a comment.

### 4.4.3 Choosing Acquisition Modes



Some acquisition modes are only available as options. Your system specific configuration might not include all modes. See  "2.5.1 Image Acquisition" P.2-11.

Before making acquisitions, you must first decide which acquisition modes to use. Programs define acquisition modes separately for Rad and Fluoro exposure. Default programs are individually preset for each performing physician by the administrator as described in  "Physicians Tab" P.17-44. Two Fluoro modes and various Rad modes exist.

| Mode   | Description   |
|--------|---|
| Fluoro | Fluoroscopy: Provides low X-ray dose images for help in positioning the catheter. |

| Mode              | Description   |
|-------------------|---|
| FluoroMap         | Fluoroscopy Map: assists in the positioning of the catheter in a complicated blood-vessel route. It highlights in white, the path of contrast agent injection. The path is retained on screen like a road map during subsequent Fluoro acquisitions to make it easy to see the moving catheter with little or no additional contrast agent. |
| DSA-MAP           | DSA-MAP: superimpose the acquired image and fluoroscopy image for observation.  |
| BlankMAP          | BlankMAP: subtract the background image and display the fluoroscopy image.  |
| Live Mode         | Live Mode: a mode not to subtract the background image.   |
| TraceMAP          | TraceMAP: superimpose the vessel contour image and fluoroscopy image for observation.   |
| DA                | Digital Angiography: the normal acquisition mode for acquiring cardiac images.  |
| DSA               | Digital Subtraction Angiography: virtually removes constant structures (flesh, bone) of no diagnostic interest, enabling enhanced blood-vessel contrast. This mode is suitable for a non-moving catheterization table, with body parts such as limbs that can be kept stationary.   |
| Flex-APS (Option) |   |
| RSM-DSA           | Real-time Smoothed Mask Digital Subtraction Angiography: (a special DSA mode), also virtually removes constant structures of no diagnostic interest, enabling enhanced blood-vessel contrast. In this mode the catheterization table can be moved and minor limb motion is acceptable.  |
| ROT-DA            | Rotational Digital Angiography: like DA but with C-arm rotation, including initial mask-acquisition rotation.   |
| ROT-DSA           | (Only for DSA option) Rotational Digital Subtraction Angiography: like DSA but with C-arm rotation, including initial mask-acquisition rotation.  |
| ROT-RSM-DSA       | Rotational Real-time Smoothed Mask Digital Subtraction Angiography: like RSM-DSA but with C-arm rotation.   |
| SPOT              | Single Image: acquisition of a single image.  |
| Stage Acquisition | Stage Acquisition: enable to acquire the image corresponding to each blood flow rate of arterial and venous phase by combining fast and slow rate.  |
| Auto Stitching    | Auto Stitching:   |

## 4 Image Acquisition

Choose a Fluoro and Rad mode by selecting a program with the desired mode in its name as follows.




You can select a program while performing a Fluoro acquisition. the selected program only applies after completion of the Fluoro acquisition.

- 1 Observe the current Rad and Fluoro program settings displayed at the bottom of the [Image Viewer] window.

Cardio / CAG[30f-10s]  
Fluo / 30pps

Default values specific to the Performing Physician are used when the new study is opened. Thereafter, the last digital program menu choice made is preserved separately for Rad and Fluoro. If the default digital program selections for Rad and Fluoro are adequate, continue with ["4.4.4 Making Live Fluoro Acquisitions" P.4-30](#).

- 2 Click Program: Fluoro] or [Program: Rad] GUI button or press the equivalent IVR NEO/IVR Shuttle/SMART Touch button.




Side Menu



IVR NEO



IVR Shuttle



SMART Touch

The appropriate digital program menu appears, either Fluoro or Rad, in a small font if launched from the GUI buttons or in a large font if launched from the IVR NEO/IVR Shuttle.

**Fluoro**  
Head ▶  
Cardio ▶  
Abd ▶  
Peri ▶  
Ablation ▶  
Other ▶  
Cancel

**Rad**  
Head ▶  
Cardio ▶  
Abd ▶  
Peri ▶  
Ablation ▶  
Other ▶  
Cancel



- 3 Once a digital program menu launches, use the keyboard, IVR NEO/SMART Touch joystick or IVR Shuttle menu selector to make program menu choices as follows.


| To move the selection highlight |                                     |
|---------------------------------|-------------------------------------|
| Keyboard:                       | Use the Up/Down arrow keys.         |
| IVR NEO/SMART Touch:            | Use the Up/Down motion of joystick. |
| IVR Shuttle:                    | Scroll Up/Down the menu selector.   |

| To change menu levels |   |
|-----------------------|---|
| Keyboard:             | Press [Enter].  |
| IVR NEO/SMART Touch:  | Use the Up/Down motion of joystick.   |
| IVR Shuttle:          | Scroll Up/Down the menu selector to highlight a menu ending with an arrow, and press the menu selector. |

| To select a menu item |                            |
|-----------------------|----------------------------|
| Keyboard:             | Press [Enter].             |
| IVR NEO/SMART Touch:  | Press the joystick button. |
| IVR Shuttle:          | Press the menu selector.   |

The highlighted menu item appears in black text with a white outline. Most parent-level menu choices lead to sub-menus as signified by menu items ending with an arrow. Menu choices that do not end in an arrow, select the indicated digital program.

For example, the Rad menu and one of its sub-menus.



You can select by directly clicking the menu with mouse instead of keyboard/joystick.

- 4 After program selection, the menu is closed. our position in the menu is remembered (Separately for Rad and Fluoro) so that the menu can be later reopened.at the same position. When you select a digital program, its name is displayed at the bottom of the [Image Viewer] window at all times except during review.



- 5 Select the parent-level to return to the parent menu and press the joystick or [Enter] key (keyboard).

Alternatively, select [Cancel] to close the menu without choosing a program.

### 4.4.4 Making Live Fluoro Acquisitions


---

Fluoroscopy provides low X-ray dose images for help in positioning the catheter. Live fluoro acquisition (no recording) can be made at any time, not just when a study is open for acquisition. To make live Fluoro acquisitions, follow this procedure.

- 1 On the Fluoro menu, select a Fluoro program.
- 2 Press the [Fluoro] switch. Live Fluoro images (Live images) are displayed.
- 3 Start injecting the contrast agent.
- 4 After acquiring the images, release the [Fluoro] switch.

If configured, the last Fluoro image acquired is held on the Acquisition monitor until another acquisition is made. And enable to perform fluoroscopy any time except opening the study for acquisition. Also, by selecting an optional image from the image selector on Acquisition monitor during fluoroscopy, enable to change the image to be displayed after fluoroscopy.

Fluoroscopy exposure time is defined by the length of time the [Fluoro] switch is pressed (with configurable maximum time). During Fluoro acquisition ([Fluoro] switch pressed), a message indicating that Fluoro acquisition is in progress is superimposed over the upper-right area of the [Image Viewer] window.



It is also possible to change the Acquisition mode from the side menu during fluoroscopy.

See  ["4.4.3 Choosing Acquisition Modes" P.4-26](#).



To change image brightness and contrast for the next Fluoro acquisition, select a fluoro loop from the Image Selector on the acquisition station and adjust image brightness and contrast using the side-menu buttons or the mouse. And then perform fluoro acquisition. The adjustments remains until a fluoro DUP is selected.



Rad (DA, DSA, RSM-DSA, ROT-DA, ROT-DSA, SCORE 3D-DA, Stage Acquisition, SPOT and SCORE CT) acquisition takes priority over Fluoro. If the [Acquisition] switch is pressed during Fluoro, the Fluoro is stopped and Rad acquisition is made. When the Acquisition switch is released, Fluoro will continue if either the Fluoro or Fluoro Record switch is still pressed.

## FluoroMAP


FluoroMAP is a mode to help positioning a catheter by superimposing blood-vessel and fluoroscopy images. Contrast agent is used to create the mask image which is subtracted from the fluoroscopy image. This mode is useful to position the catheter in a complicated blood-vessel route during study.

To make FluoroMap acquisitions, follow this procedure.


1 Select DUP for any purpose from Fluoro program.

2 Click [MAP]-[Fluoro] on the side menu.

A message of current mode (FluoroMAP MASK creation) is displayed at the bottom-left of the acquisition monitor.



Acquisition Monitor




Reference Monitor

3 Press [Fluoro] switch. And start acquiring the background image.

4 Mask creation will start after acquiring the background image. And start injecting contrast agent.

5 When the mask image for catheter operation is acquired, release [Fluoro] switch.

A message of current mode (FluoroMAP) is displayed at the bottom-left of the acquisition monitor.



6 When keep pressing [Fluoro] switch, the subtraction image of mask and Fluoro image created in step 4 is displayed.

7 FluoroMAP images can be recorded.

👉 ["4.4.5 Recording Fluoro" P.4-38.](#)

- 8 Select [MAP]-[Vessel Power] to change mask image brightness from ACQ PC side menu. Right-click the [Vessel Power] button to increase and left-click to decrease the image brightness.




- 9 If the desired effect is not achieved, repeat step 2 to create the mask again.
- 10 To return to normal fluoro mode, select [MAP]-[ON] (MAP mode) on the side menu. Keep FluoroMAP mask until new mask is set.
- 11 During normal (non-Map) fluoro mode, click [MAP]-[OFF] (MAP mode) on the side menu to change the mode to MAP mode.

### ■ DSA-MAP

DSA-MAP is a mode to subtract the background image and observe by superimposing existing blood-vessel and fluoroscopy images while positioning the catheter. As it is not necessary to acquire the blood-vessel image again, radiation dosage and volume of contrast agent is reduced.

To make DSA-MAP acquisitions, follow this procedure.


- 1 Select a DUP for any purpose from Fluoro menu.
- 2 Select the blood-vessel image from the Image Selector on the Reference or Acquisition monitor.  
Cannot use a vessel contour image which is saved as a still image.
- 3 If entire blood-vessel image is not necessary, click the [Region] button on the DSA Tools bar and set ROI for the necessary blood-vessel.



- 4 Click [MAP]-[DSA] from Reference monitor side menu and set the selected image in Step 2 as the mask image. To set the image of Acquisition monitor, click [MAP]-[DSA] on the Acquisition monitor.


A message of current mode (DSA-MAP) is displayed at the bottom-left of the Acquisition monitor.



- 5 When pressing the [Fluoro] switch, start acquiring the background image.
- 6 Superimposed image of Mask image created in Step 4 and the fluoroscopy image is displayed.
- 7 DSA-MAP images can be recorded.  
 ["4.4.5 Recording Fluoro" P.4-38.](#)
- 8 [MAP]-[Vessel Power] to change mask image brightness from ACQ PC side menu. Right-click the [Vessel Power] button to increase and left-click to decrease the image brightness.



- 9 To return to normal mode, select [MAP]-[ON] (MAP mode) on the side menu. Keep DSA-MAP mask until new mask is set.
- 10 During normal (non-Map) Fluoro mode, click [MAP]-[OFF] (MAP mode) on the side menu to change the mode to MAP mode.





### BlankMAP


BlankMAP is a mode to subtract the background image such as bone and display fluoroscopy images. This is useful when checking device and medical agent added after BlankMAP is started.

To make BlankMAP acquisitions, follow this procedure.

- 1 Click [MAP]-[Blank] from side menu on the reference monitor. A message of current mode (BlankMAP) is displayed at the bottom-left of the acquisition monitor.




- 2 Press [Fluoro] switch. And start acquiring the background image.
- 3 Acquisition of subtraction image will start after acquiring the background image.
- 4 BlankMAP images can be recorded.  
 ["4.4.5 Recording Fluoro" P.4-38.](#)
- 5 To acquire the background image again, click [MAP]-[Reset] and press [Fluoro] switch.



Reference Monitor

- 6 To return to normal fluoro mode, click [MAP]-[ON] (MAP mode) on the side menu.




Reference Monitor

## LIVE Mode

LIVE mode is a mode not to subtract the background image such as bone. This mod is available in FluoroMAP and DSA-MAP acquisition.


- 1 Click [MAP]-[OFF] (MAP mode) on side menu of the reference monitor. If MAP mode is not valid, perform FluoroMAP or DSA-MAP.

 ["FluoroMAP" P.4-31](#) and ["DSA-MAP" P.4-32](#)


- 2 Click [MAP]-[LIVE] from the side menu. A message of current mode (Live) is displayed on the bottom-left of the acquisition monitor.




Acquisition Monitor



Reference Monitor


- 3 Press [Fluoro] switch. Acquire fluoroscopy image without subtracting the background image.
- 4 Live mode images can be recorded.  
 ["4.4.5 Recording Fluoro" P.4-38](#)
- 5 Click [MAP]-[Live] on the side menu to return to MAP mode.



Reference Monitor



If [Coordination with Geometry] is available on [Physicians] tab, normal fluoroscopy may be performed though it is set to MAP mode.

 [Selecting Geometry Configuration] of "[Physicians Tab](#)" P.17-44.

- Changing C-arm Angle (RAO/LAO/Cranial/Caudal) and SID:  
When changing the C-arm angle and SID, normal fluoroscopy is performed though it is set to MAP mode. And when return to the original position, MAP will be performed.
- Changing the C-arm Parallel Position/Bed Position:  
When changing the position of C-arm or bed, mask image will be deleted and MAP mode will be canceled. It is set to normal fluoroscopy mode, and MAP will not be performed though return to the original position.



- In case of LIDM, a fluoroscopy image during radiation and vascular MAP image can be displayed at the same time.
- Vascular MAP image is displayed on the Acquisition monitor during fluoroscopy and fluoroscopy image is displayed on Reference monitor in Examination room.
- When adding LIVE monitor, fluoroscopy image can be displayed on the LIVE monitor.
- When releasing the [Fluoro] switch, displayed LIVE image will be switched to currently selected image on Reference and Acquisition monitor. At this point, it requires a little bit of time for switching. And nothing will be displayed on the LIVE monitor except for performing fluoroscopy.
- Image quality of LIVE image is a little bit different from a normal fluoroscopy image.
- LIVE image flips with HV flip on the Acquisition monitor during fluoroscopy.
- Zoom cannot apply to the LIVE image.



## TraceMAP

TraceMAP is a mode that fluoroscopy image and vessel contour image generated with DSA mode can be displayed as a superimposed. Interior of a blood-vessel is displayed in a transparent mode, so this is available when treat a large vessel.

To make TraceMAP acquisitions, follow this procedure.

- 1 Select DSA radiography program.
- 2 Press the [Radiography] switch and inject contrast agent to acquire a blood-vessel image.
- 3 After the acquisition, release the switch.
- 4 Select the blood-vessel image from the Image Selector on the Reference monitor. Or use the Peak Hold function to display the whole blood-vessel image.

 ["9.4 Peak Hold" P.9-6](#)

- 5 Select the image to enhance contour and click the [Contour] button on [DSA Tools] bar.
- 6 Apply contour enhancement with default parameters for the blood-vessel image on the monitor. Adjust parameters on the dialog and display an appropriate image.

 ["9.5 Contour Enhancement" P.9-8](#)




- 7 Enable to replay the back ground image while displaying a contour image. Confirm that the vessel is enhanced correctly.
- 8 If the entire contour enhancement image is not necessary, click [Region] button on the [DSA Tools] bar and set ROI for the necessary contour.


 ["9.6 Mask Region" P.9-11](#)

- 9 Click [MAP]-[DSA] on the side menu and select created contour enhancement image as a superimposed image. If contour is displayed, the DSA-MAP will operate as the [TraceMAP] button.

A message of current mode (TraceMAP) is displayed at the bottom-left of the Acquisition monitor.



Acquisition Monitor



Reference Monitor

- 10 Press the [Fluoro] switch to display superimposed image of blood-vessel contour and fluoroscopy images. Refer to step 5 to 9 of making DSA-MAP Acquisitions for subsequent operations.

### 4.4.5 Recording Fluoro

Although the main purpose of the system is to record Rad acquisitions, it is also possible to record Fluoro acquisitions in one of two modes: Direct Record or Last N Seconds. Direct Record requires you to press a button for each Fluoro recording. Last N Seconds causes the most recent Fluoro images to always be recorded in a temporary buffer. The images in the buffer can then be saved when desired. The Fluoro record mode is set individually for the default physician and each performing physician.

 ["Physicians Tab" P.17-44](#)



The images displayed after fluoroscopy can be changed by selecting any loop from the image selector on the Acquisition monitor during fluoroscopy.

 ["The acquisition procedures are describe in the following sections." P.4-24](#)

#### Direct Fluoro Record




To use Direct Record, the default physician must be configured as [Fluoro Record]-[Direct] or if you are logged in as a configured performing physician, you must be configured the same.

 ["Physicians Tab" P.17-44](#)

Unable to record the LIVE fluoroscopy image if it is displayed with the optional feature of LIDM.

Unable to use Direct Record during fluoroscopy without a record (while pressing a [Fluoro] switch).

To use Direct Fluoro Record, follow this procedure.

- 1 Choose your Rad and Fluoro programs as described in  "4.4.3 Choosing Acquisition Modes" P.4-26.
- 2 Make Fluoro acquisitions.
- 3 When it is desired to record the next Fluoro acquisition, click [Fluoro]-[Direct Record] button on the GUI side menu.

The GUI button appears pushed in. Press the [Fluoro] switch to begin Fluoro acquisition with recording. Release the [Fluoro] switch to stop acquisition and recording.

The GUI [Fluoro]-[Direct Record] button pops back out.

4


#### Last N Seconds Fluoro Record



To use Last N Seconds, the default physician must be configured as [Fluoro Record]-[Last N Seconds] or if you are logged in as a configured performing physician, you must be configured the same.

 "Physicians Tab" P.17-44

To use Last N Seconds Fluoro Record, follow this procedure.

- 1 Choose your Rad and Fluoro programs as described earlier in  "4.4.3 Choosing Acquisition Modes" P.4-26.
- 2 Make Fluoro acquisitions.
- 3 If it is desired to save your most recent Fluoro acquisition, click [Fluoro]-[Save] button on the GUI side menu or press the IVR NEO/IVR Shuttle/SMART Touch [Save Fluoro] button.

The Fluoro acquisition is saved and transferred to the archive server and Reference computer. Optionally, click the [Fluoro]-[View] button on the GUI side menu to first playback the recorded (but not yet saved) Fluoro images on the Acquisition monitor to see if they are worth saving.




If the Fluoro acquisition duration is longer than the configured maximum number of seconds (shown in program name), only the first portion (up to configured limit) is kept in the recording buffer and made available for saving.

### 4.4.6 Recording Rad

---

Digital Angiography is the normal acquisition mode for acquiring cardiac images. Rad acquisitions are always recorded.

Rad exposure time is defined by the length of time the [Acquisition] switch is pressed (with configurable maximum number of images). During Rad acquisition ([Acquisition] switch pressed), a message indicating the active mode is superimposed over the upper-right area of the Image Viewer.



Playback of the acquired loop starts.



Acquisition (DA, DSA, RSM-DSA, ROT-DA, ROT-DSA, SCORE 3D-DA, SCORE 3D-DSA, Stage Acquisition, SPOT and SCORE CT) takes priority over Fluoro. If [Acquisition] switch is pressed during Fluoro, the Fluoro is stopped and DA acquisition is made. When the [Acquisition] switch is released, Fluoro will continue if either the [Fluoro] or [Fluoro Record] switch is still pressed.



If required, choose a comment as described in ["4.4.2 Adding a Comment" P.4-25](#).

To make acquisitions, follow this procedure.

#### DA

DA is the normal acquisition mode for cardiac images.

- 1** Select a DA program on the Rad menu.
- 2** Press the [Acquisition] switch.
- 3** X-ray exposure starts and live images are displayed and recorded. Start injecting the contrast agent.
- 4** When image acquisition is complete, release the [Acquisition] switch.
- 5** Playback of the acquired loop starts.

## DSA

*If the catheterization table must be moved, or the limb cannot be kept motionless, use RSM-DSA instead (described in the next section).*

DSA mode virtually removes constant structures (flesh, bone) of no diagnostic interest, enabling enhanced blood-vessel contrast. These modes are suitable for a non-moving catheterization table, with motionless body parts such as limbs. DSA acquisitions are always recorded.

- 1 Select a DSA program on the Rad menu.
- 2 Press the [Acquisition] switch.
- 3 X-ray exposure starts and the first images are used to create the mask, as configured in the DUP.
- 4 Start injecting the contrast agent for subtraction. Images are displayed subtracted to reveal enhanced blood-vessel contrast.
- 5 When image acquisition is complete, release the [Acquisition] switch.
- 6 Playback of the acquired loop starts.

## Flex-APS (Option)

Flex-APS is a special DSA mode virtually adjusts the position of image automatically to reduce the effect of body movement when creating a subtraction image.



Do not use Flex-APS for a region with movement inside the body such as flatus. Artifact may be occurred.

- 1 Select a Flex-APS program from Rad program.
- 2 Press [Acquisition] switch.
- 3 X-ray exposure starts and the mask image is displayed. And then start injecting the contrast agent for subtraction.
- 4 When image acquisition is complete, release the [Acquisition] switch.
- 5 Playback of the acquired loop starts.



Enable to OFF the acquired automatic position correction for acquired image.

Refer to  ["9.3 Reregistration" P.9-5](#) for details.

### RSM-DSA

RSM-DSA, a special DSA mode, also virtually removes constant structures (flesh, bone) of no diagnostic interest, enabling enhanced blood-vessel contrast. In this mode, it is acceptable to have catheterization table motion and/or minor limb motion. RSM-DSA acquisitions are always recorded.

- 1 Select RSM-DSA program on the Rad menu.
- 2 Press the [Acquisition] switch.
- 3 Start injecting the contrast agent. Images are displayed subtracted to reveal enhanced blood-vessel contrast. It is possible to move the catheterization table during image acquisition.
- 4 When the image acquisition is complete, release the [Acquisition] switch.
- 5 Playback of the acquired loop starts

### ROT-DA and ROT-RSM-DSA

Acquisitions can be made by rotating C-arm.

- 1 In case of Bi-plane mode, move the Lateral C-arm to the park position.
- 2 Select a DA or RSM-DSA program on the Rad menu.
- 3 A message is displayed on the top of the Acquisition monitor.
- 4 Register the C-arm rotation start and stop position as follows:
  - 1 Position the C-arm to the stop position first and press the left-arrow button on the C-arm local console (or remote console).
  - 2 Position the C-arm to the rotation start position and press the blinking button again to register the position.
- 5 While the C-arm is close to the catheterization table, a message is displayed on the Acquisition monitor to press the [SET] button on the C-arm local console (or remote console.)

Press the [SET] button on the C-arm local console (or remote console) to position the C-arm where the registration is enabled.

- 6 After registration, press the [TEST RUN] button on the C-arm local console (or remote console) to confirm the position.



The [TEST RUN] uses fluoroscopy to confirm that the focus region stays the center of the field of view while rotating. The test run ends when the button is released.

- 7 To position the C-arm after test run completion, press the [CANCEL] button to cancel the rotation mode.
- 8 Press and hold the [SET] button to move the C-arm to the start position. Hold the [SET] button until the indicator of the switch stops blinking and illuminates.  
When the system is ready for acquisition, the message "X-RAY READY" is displayed on the Acquisition monitor.
- 9 Prepare the injector, and press and hold the [Acquisition] switch to perform radiography.
- 10 The acquisition ends when the [Acquisition] switch is released or when the maximum number of acquired images is reached (acquisition continues even after the C-arm reaches the rotation stop position).  
During acquisition, a message indicating the active mode is superimposed over the upper-right area of the Image Viewer window.

**ROT DA**

**ROT RSM**

#### **ROT-DSA**

- 1 In case of Bi-plane mode, move the Lateral C-arm to the park position.
- 2 Select a ROT-DSA program on the Rad menu.
- 3 A message is displayed on the top of the Acquisition monitor.
- 4 Register the C-arm rotation start and stop position as follows:
  - 1 Position the C-arm to the stop position first and press the left-arrow button on the C-arm local console (or remote console).
  - 2 Position the C-arm to the rotation start position and press the blinking button again to register the position.

- 5 While the C-arm is close to the catheterization table, a message is displayed on the Acquisition monitor to press the [SET] button on the C-arm local console (or remote console.)

Press the [SET] button on the C-arm local console (or remote console) to position the C-arm where the registration is enabled.

- 6 After registration, press the [TEST RUN] button on the C-arm local console (or remote console) to confirm the position.



The [TEST RUN] uses fluoroscopy to confirm that the focus region stays the center of the field of view while rotating. The test run ends when the button is released.

- 7 To position the C-arm after test run completion, press the [CANCEL] button to cancel the rotation mode.
- 8 Press and hold the [SET] button to move the C-arm to the start position. Hold the [SET] button until the indicator of the switch stops blinking and illuminates.  
When the system is ready for acquisition, the message "X-RAY READY" is displayed on the Acquisition monitor.
- 9 Start injecting the contrast agent and press the [Acquisition] switch. The C-arm starts rotating and mask acquisition starts.
- 10 Acquisition ends when the [Acquisition] switch is released, when the C-arm reaches the rotation stop position or when the maximum number of acquired images is reached.
- 11 After acquiring the mask images, press the [SET] button to move the C-arm back to its start position.
- 12 Press the [Fluoro] switch.  
The MAP image is displayed on the Acquisition monitor.
- 13 Ask a patient to hold his/her breath at the position of the least misregistration. Prepare the injector, and press and hold the [Acquisition] switch to perform radiography.
- 14 The acquisition ends when the C-arm reaches the position in which the mask acquisition has ended or the [Acquisition] switch is released.





If DSA image acquisition occurs for a smaller rotational range than mask acquisition, then any mask images that do not have a corresponding live image are not used.

## SPOT

SPOT acquisition is the normal acquisition mode for acquiring cardiac images. DA (Rad) acquisitions are always recorded. To make SPOT acquisitions, follow this procedure.

- 1 Select a SPOT program on the Rad menu.
- 2 If required, choose a comment as described in ["4.4.2 Adding a Comment" P.4-25](#).
- 3 Press the [Acquisition] switch. X-ray exposure starts and one live image is displayed and recorded.
- 4 Release the [Acquisition] switch.

During SPOT acquisition ([Acquisition] switch pressed), a message indicating the active mode is superimposed over the upper-right area of the Image Viewer window.



SPOT acquisition takes priority over Fluoro. If the [Acquisition] switch is pressed during Fluoro, the Fluoro is stopped and SPOT acquisition is made. When the [Acquisition] switch is released, Fluoro will continue if either the [Fluoro] or [Fluoro Record] switch is still pressed.

## Precession/pendulum DA and RSM-DSA

To make Precession/Pendulum DA or RSM-DSA acquisitions, follow this procedure.

- 1 Select a DA or RSM-DSA program on the Rad menu.
- 2 If required, choose a comment as described in ["4.4.2 Adding a Comment" P.4-25](#).
- 3 On the C-arm local console (or remote console), press the [Precession] or [Pendulum] switch.
- 4 On the C-arm local console (or remote console), press the [SET] button until the C-arm has fully moved to the default (start) position.
- 5 Perform fluoro acquisitions as needed.

- 6 When ready, press the [Acquisition] switch.  
The C-arm motion and acquisition starts.
- 7 Inject contrast agent as needed.
- 8 Acquisition continues until the [Acquisition] switch is released or the C-arm moves a total of three cycles.  
When acquisition is complete, playback of the acquired loop starts.  
During acquisition, a message indicating the active mode is superimposed over the upper-right area of the Image Viewer window.



### ■ Stage Acquisition

On this mode, acquire DSA image by reducing acquisition rate for each specified number of frames that are set within the radiography program, in one acquisition. Enable to acquire the image corresponding to each blood flow rate of arterial and venous phase by combining fast and slow rate.

- 1 Select the stage acquisition mod program from radiography program.
- 2 Press the radiography switch.
- 3 X-ray is exposed and the mask image is displayed. Subtraction is started, so inject the contrast medium.
- 4 Acquisition rate is reduced after each scheduled number of frames are acquired. Acquisition will be completed when specified number of images are acquired.
- 5 Automatic review of acquisition data is started after radiography.

### 4.4.7 Flipping Images

Images can be flipped horizontally or vertically at acquisition time. Click the [H Flip] button (left in the following figure) or [V Flip] button (right in the following figure) on the side menu of Acquisition monitor to flip the images for the next acquisition horizontally or vertically, respectively.



Default setting is determined depends on a patient's orientation set at the start of study.

| Orientation | H Flip | V Flip | Orientation | H Flip | V Flip |
|-------------|--------|--------|-------------|--------|--------|
| HFS         | OFF    | OFF    | FFS         | ON     | ON     |
| HFP         | ON     | OFF    | FFP         | OFF    | ON     |
| HFDR        | OFF    | OFF    | FFDL        | OFF    | ON     |
| HFDL        | OFF    | OFF    | FFDR        | OFF    | ON     |

You can flip fluoroscopy or radiography images, but only for the next acquisition.

To flip the acquired image, click the [Filter]-[H Flip] or [Filter]-[V Flip] button on the side menu of Reference monitor.

["Applying Filters" P.7-7](#)


Flip fluoroscopy Live image according to the HV flip of Acquisition monitor during fluoroscopy.

### 4.4.8 Changing the Acquisition Parameter Settings


The injection delay settings allows for the contrast agent to spread the area under study and include the injection-acquisition delay and the mask-acquisition delay in advance. And enable to change values during study. The injection-acquisition delay is the delay in seconds between the activation of the acquisition switch and the beginning of acquisition. The mask-acquisition delay is the delay in seconds between the end of mask acquisition and the beginning of live acquisition.

Default values are configured for each Fluoro and Rad program. To temporary change these values during study, follow this procedure.

- 1 Click [Acq param] on the side menu of Acquisition monitor.



Single-plane



Bi-plane

The [Set Acquisition Parameters] window opens. The default values are displayed.

[If DUP other than HQ-DSA is selected]

The 'Set Acquisition Parameters' window shows the 'Acquisition' section with '15 fps' and 'Total: 151 frames (10 sec)'. The 'Rad Time' is set to '10'. The 'Delay' section shows 'Rad I-A (s): 0' and 'Rad M-A (s): 0'. 'Apply' and 'Close' buttons are at the bottom right.

The 'Set Acquisition Parameters' window shows the 'Acquisition' section with '15 fps' and 'Total: 234 frames (15 sec)'. The 'Rad Time' is set to '15' and 'Rad Mask' is set to '8'. The 'Delay' section shows 'Rad I-A (s): 0' and 'Rad M-A (s): 0'. 'Apply' and 'Close' buttons are at the bottom right.

[If DUP of HQ-DSA is selected]

The 'Set Acquisition Parameters' window shows the 'Stage Acquisition' section with '6 fps', '4 fps', '3 fps', '2 fps', '1 fps', 'Rad Mask', and 'Total: 56 frames (30 sec)'. The 'Rad Time' is set to '3', '0', '2', '2', '23', and '1'. The 'Delay' section shows 'Rad I-A (s): 0' and 'Rad M-A (s): 0'. 'Apply' and 'Close' buttons are at the bottom right.

- 2 Change the parameters as needed and click [Apply] button.

Delay units are in seconds and can be set to a precision of up to one decimal. These values remain until you select a Fluoro or Rad program.

["4.4.3 Choosing Acquisition Modes" P.4-26](#)

- 3 Click the [Close] button to close the window.

- 4 You are now ready to begin acquisition.



Injection-Acquisition Delay can be displayed at the bottom-left of the Acquisition and Reference monitors.

["Reference Monitor" P.17-17](#)















#### 4.4.9 Manually Resetting the Dose

---

Automatic dose reset can be configured to occur at the end of a study or the beginning of the next one. If a study is not open, manual dose reset is also available. To manually reset the dose, click [Funct]-[Reset] on the side menu of Acquisition monitor.

## 4.5 Actions When Study Open for Acquisition

Although the purpose of having a study open for acquisition is to acquire loops, other actions can be performed.

| No. | Title   | Reference  |
|-----|---|--|
| 1   | Continuing a Study                                      |  <a href="#">"4.5.1 Continuing a Study" P.4-50</a>                                      |
| 2   | Using Select  |  <a href="#">"4.5.2 Using Select" P.4-51</a>  |
| 3   | Using Split   |  <a href="#">"4.5.4 Filtering Images by C-arm Position" P.4-53</a>                      |
| 4   | Filtering Images by C-arm Position                      |  <a href="#">"4.5.4 Filtering Images by C-arm Position" P.4-53</a>                      |
| 5   | Moving C-arm to Match Selected Image (Reposition)       |  <a href="#">"4.5.5 Moving C-arm to Match Selected Image (Reposition)" P.4-54</a>       |
| 6   | Moving Table to Match Selected Image (Table Reposition) |  <a href="#">"4.5.6 Moving Table to Match Selected Image (Table Reposition)" P.4-55</a> |
| 7   | Sending Selected Image                                  |  <a href="#">"4.5.7 Sending Selected Image" P.4-55</a>                                  |
| 8   | Deleting Selected Image                                 |  <a href="#">"4.5.8 Deleting Selected Image" P.4-56</a>                                 |
| 9   | Displaying Previous Studies                             |  <a href="#">"4.5.9 Displaying Previous Studies" P.4-56</a>                             |
| 10  | Sketch Display  |  <a href="#">"4.5.10 Sketch Display" P.4-57</a>  |
| 11  | Working with Reference Images                           |  <a href="#">"4.6 Working With Reference Images" P.4-60</a>                           |
| 12  | Annotating  |  <a href="#">"7.7 Working with Image Annotations" P.7-17</a>                          |
| 13  | Saving  |  <a href="#">"7.8 Saving Images" P.7-20</a>   |
| 14  | Printing  |  <a href="#">"7.9 Print Images" P.7-21</a>  |
| 15  | Analyzing   |  <a href="#">"8 Performing Angiographic Analysis" P.8-1</a>                           |
| 16  | Adjusting DSA Images                                    |  <a href="#">"9 DSA Image Adjustment Tools" P.9-1</a>                                 |

### 4.5.1 Continuing a Study

Studies initiated from a worklist server can be continued with additional acquisitions even after closure. To continue such a study, click [New] on the side menu and then find and select the desired study in the worklist. Perform additional acquisitions and create additional images as desired. The new loops and images are added at the end of the study.

### 4.5.2 Using Select


The [Select] function is provided primarily for use by the physician in the examination room to enable easy loop selection with the IVR NEO/SMART Touch buttons and joystick. Additionally, the control-room technician may use equivalent mouse actions or keyboard keys.

Use [Select] as follows.

- 1 Press the [16-Up Selector] button on the IVR NEO/SMART Touch or the [Funct]-[Select] button on the side menu.

The Image Viewer changes to a 4x4 grid of image thumbnails. The loop number and C-arm position CRA-CAU/RAO/LAO) are superimposed in large letters over each thumbnail. Also, a Loop/Still Image icon is superimposed over the upper-right corner of each thumbnail, indicating whether a loop or still image is represented.

A light-colored selector box surrounds the thumbnail of the current loop, the loop that was on screen when [Select] was chosen.



- 2 To choose a particular loop to display at full size, first move the selector (light-color box) with the IVR NEO/SMART Touch joystick, keyboard arrow or [PageDown]/[PageUp] keys, so that the thumbnail of the desired loop is surrounded.
- 3 Click the IVR NEO/SMART Touch joystick button or press [Enter] on the keyboard. Alternatively, just click any thumbnail with the mouse.

The selected loop appears full size in the Image Viewer window.

- 4 If there are more than 16 loops and still images in the study, scroll the thumbnail list downward by moving the IVR NEO/SMART Touch joystick down or by right-click upward dragging the Additional thumbnail grid with the mouse, or by pressing the keyboard down arrow or [PageDown] key. Additional thumbnails are revealed from below and the initial thumbnails are scrolled off the top.

Select a loop for display as already described.

- 5 To exit the Select mode without choosing a loop for display, press [16-Up Selector] button on the IVR NEO/SMART Touch or the [Funct]-[Select] button on the side menu.

### 4.5.3 Using Split

(For Reference monitor only.)

To split the Image Viewer into two side-by-side panes, press the [Single Display]/[Split Display] button on the IVR NEO/SMART Touch or click [Funct]-[Split] on the side menu. The Image Viewer splits down middle with the selected loop displayed on the left pane and the right pane empty. The next loop will be displayed in the right pane. A selection box surrounds the left pane.

To display a loop in the right pane, press the [Left/Right Focus] button on the IVR NEO/SMART Touch or click the right half of the Image Viewer window (the selection box now surrounds the right pane) and then choose a loop (or image) in the Image Selector.

When both panes contain loops, the loop in the selected pane can be played and the loop in the other pane is paused. Pan to the left or right to see the entire loop.




Image processing (Brightness, Contrast, Auto Window Level, Filters, LUT) and playback controls are only applicable to the selected pane.




To evaluate different image processing settings, place the same loop in both panes, select a pane and make adjustments. The effects are shown only in the selected pane. This way, you can experiment with different image processing to see how the loop looks both with and without the image processing.

To un-split the display, choose the Split feature again. Whatever was in the last-selected pane is shown in the full Image Viewer window.

#### 4.5.4 Filtering Images by C-arm Position


To filter the Image Selector so only icons for images that were acquired with a C-arm position that matches the current C-arm position (within configurable tolerances), follow this procedure.

- 1 Position the C-arm at the desired position.
- 2 Press the IVR NEO/SMART Touch [Filter by C-Arm] button. (LED of IVR NEO turns on)
- 3 Alternatively, click the Image Filtering button on the bottom GUI bar and then, in the Image Filtering list, choose [Match C-Arm].




The Image Selector is filtered so only the icons of images that were acquired at the current C-arm position are shown with the LAO/RAO and CRA/CAU angles at the top like this.

The Image Selector is no-longer filtered; icons are shown like this.




- 4 To revert an unfiltered Image Selector, press the IVR NEO/SMART Touch [Filter by C-Arm] button again.

LED of IVR NEO button turns off, or click the Image Filtering button on the bottom GUI bar and then, in the Image Filtering list, choose [All Images].




The Image Selector is no-longer filtered; icons are shown like this.



### 4.5.5 Moving C-arm to Match Selected Image (Reposition)

To move the C-arm to the position at which the selected image was acquired, follow this procedure.

- 1 Select the desired image or loop.
- 2 Press the IVR NEO/SMART Touch [C-arm Reposition] button (LED of IVR NEO turns on), or click [Funct]-[Move] on the side menu.



- 3 Press the C-arm [Direct Memory] button as it blinks.

The C-arm moves to the position at which the selected image was acquired. If the C-arm cannot move to the requested position, a message to that effect is displayed.

#### 4.5.6 Moving Table to Match Selected Image (Table Reposition)




Enable to use when angiographic catheterization table KS-100 is installed.

4

To move the table to the position at which the selected image was acquired, follow this procedure.

- 1 Select the desired image or loop.
- 2 Press the IVR NEO/SMART Touch [Table Reposition] button (LED of IVR NEO turns on), or click [Func]-[Move] on the side menu.




- 3 Press the table [SET] button as it blinks.

The table moves to the position at which the selected image was acquired. If the table cannot move to the requested position, a message to that effect is displayed.

#### 4.5.7 Sending Selected Image

Enable to send the selected image. To send the image, follow this procedure.

- 1 From image selector on Reference monitor, put a check mark in the check box at the bottom left of a desired still image or loop.  
Right-click on the image selector and select [Select All] on the displayed pop-up menu, enable to select all the images.

- 2 Click [Func]-[Send] button on the side menu. A pop-up menu to select destination is displayed. Or right-click the image selector and click [Send], a pop-up menu to select destination is displayed.
- 3 Check [Process Before Transfer] in a pop-up menu for image processing.  
See  "DICOM Tab" P.17-24 for processing setting.
- 4 Click the destination.  
Replaying frame will be sent.

### 4.5.8 Deleting Selected Image

---

Enable to delete the selected image. To delete the image, follow this procedure.

- 1 From image selector on Reference monitor, put a check mark in the check box at the bottom left of a desired still image or loop.  
Right-click on the image selector and select [Select All] on the displayed pop-up menu, enable to select all the image.
- 2 Right-click on the image selector and select [Delete] on the pop-up menu, and the message is displayed.


### 4.5.9 Displaying Previous Studies

---

While a study is open for acquisition, the performing physician can request that the control room technician displays an earlier study for the same patient (same Patient ID) on the Reference monitor. The technician can then find the specific image or loop that the physician wants to see for comparison purposes.

To load a previous study for the patient, follow this procedure.

- 1 Click [Open] on the side menu.



The Studies Management window appears and a search is performed for all studies with the same Patient ID as the active study.


- 2 Double-click the desired study or select it and click [View]. You are prompted with the following message:  
*"Do you want the study STUDY NAME\_ACCESSION NUMBER to a replace the current study?"*
- 3 Click [Yes] to replace whatever is displayed with the selected study.



This only refers to temporarily replacing the images displayed on the Reference monitor with the selected study. Later, you can re-synchronize on the active study in acquisition by clicking [Sync] on the side menu.

The study is displayed in the Image Viewer. The physician in the examination room can now compare the current images against the previous. The current images are shown on the Acquisition monitor, and the previous on the Reference monitor.

- 4 Loops can be played back and specific frames can be shown on the both the Acquisition and Reference monitors via keyboard keys as described in ["5.2.1 Controlling Image Playback" P.5-5](#).
- 5 Once image comparison is finished, click [Sync] on the side menu for re-synchronize the Reference monitor with the active study on the Acquisition monitor.





If a new acquisition is not taking place, and the Acquisition monitor is playing the last acquired loop, the Reference monitor will begin playing the same loop (if configured).

#### 4.5.10 Sketch Display

Use the Sketch Display on Acquisition monitor.

Enable to place a sketch for such as a stent. To display the sketch on the image, follow this procedure.

- 1 On Acquisition monitor, display a frame such as DSA image, which you want to place a sketch.
- 2 Click [Funct]-[Sketch Ed.] on the GUI side menu. [Edit] window is displayed. Straight line is selected as a default.





### NOTE

Enable to set a type of default line selected on [Edit] window.

See [Selecting 2D Guidance Configuration] on ["Physicians Tab" P.17-44](#).

- 3 Click on the image. A point is displayed. And click another place to display a straight line which links these points. The straight line is specified in pair and the third point is the starting point of new straight line.
- 4 Polyline, Curve, and Free Line can be placed also. In this case, after selecting a desired sketch from Sketch Edition window, click on the image.
- 5 After placing the sketch, close the Sketch Edition window. The sketch continues to display after closing the window.
- 6 Click [Funct]-[Sketch] on the side menu of the Acquisition monitor to hide the sketch, or follow the instruction below.
  - If the mode is MAP mode, click [MAP]-[ON] (setting MAP mode) on the side menu and set to fluoroscopy mode.
  - On [Physicians] tab, change the geometry when [Coordination with Geometry] is set.



### NOTE

The sketch switches Display and Hide as follows depending on the changed geometry.

- C-arm Angle (RAO/LAO/Cranial/Caudal) and SID.

If C-arm angle and SID are changed, the sketch will be hidden, but displayed again if return to original position. See [Selecting 2D Guidance Configuration] on ["Physicians Tab" P.17-44](#).


- C-arm Parallel Position and Catheterization Table Position

If C-arm or catheterization table position is changed, the sketch will be deleted at this point and not be displayed if return to original position.

### 4.5.11 Circle Display


Use the circle to align the subject with X-ray irradiation area and the center of the image.

- 1 From the side menu in the Acquisition monitor, select [Funct.]-[Circle].



For bi-plane system, select [Circle] button on both frontal and lateral Acquisition monitor.

- 2 Make sure the circles appear on the image in Acquisition monitor.



- 3 Align the subject to the center of the circle.

## 4.6 Working With Reference Images



Reference images can be added or deleted for the study that is open for acquisition. Later, in review, reference images can be viewed as still images.

### 4.6.1 Introduction

Any image loop frame can be selected as the reference image for the loop. At study closure, the reference images are created and saved with the study as still images.


### 4.6.2 Adding a Reference Image

To add a reference image to a study that is open for acquisition, follow this procedure.

- 1 Use the IVR NEO/IVR Shuttle/SMART Touch [Monitor Select] button to choose the Reference or Acquisition monitor.  
Move the joystick up or down to jump between loops and then pause on the image you wish to select as a reference image for the loop. If you wish to choose a loop that does not yet have a Reference image, make sure that the IVR NEO is not in Reference mode (IVR NEO [Reference Mode] button is OFF.) If configured, the indicator "Ref" also appears in the upper-left corner of each monitor that is Reference mode.
- 2 Alternatively, click the desired Image Selector icon and pause on the image you wish to select as a reference image for the loop.
- 3 Press the IVR NEO/IVR Shuttle [Add Reference Image] button or SMART Touch [Add Ref Image], or on the keyboard, press the [Insert] key.

The selected image frame is marked for saving in the study as a full-size reference image when the study is closed. The Image Selector icon is replaced with a small version of the selected image.

As a reminder, the image frame number of the selected image is added inside square brackets to the right of the loop number like this (image frame 39 of loop 12).



The reference images themselves, do not appear in the Image Selector until the study is closed and opened in review. When reviewing, scroll down in the Image Selector to see the reference images.





All reference image selections (image frame number in square brackets on Image Selector icon) on the Reference monitor are automatically shown on the Acquisition monitor.

### 4.6.3 Displaying Images Selected as Reference

Since the actual reference images are only created at study closure, they are not yet available for display. However, you can jump to each loop that has a selected reference image and automatically pause on the selected image.

With the IVR NEO/SMART Touch, select the desired monitor, either Reference or Acquisition. Ensure that the IVR NEO/SMART Touch is in Reference mode, pressing the IVR NEO/SMART Touch [Reference Mode] button if it is not already ON. Move the joystick up/down or use the [Select Reference Image] button on the IVR Shuttle to jump to the previous/next loop with a reference image and automatically pause on the selected image.



The Reference and Acquisition monitors can be independently set to Reference mode. If configured, the indicator "Ref" also appears in the upper-left corner of each monitor that is in Reference mode. The indicator disappears when the monitor is not in Reference mode.

Alternatively, regardless of Reference mode state, press the [Page Up/Page Down] keyboard keys to jump to the previous/next loop with a reference image.

### 4.6.4 Deleting a Reference Image

Move the IVR NEO/SMART Touch joystick up/down or use the [Select Reference Image] button on the IVR Shuttle to select the desired loop.

Alternatively, on the keyboard, press the [Page Up/Page Down] keys to jump to select the desired loop and press the keyboard [Delete] key.


The image frame number in square brackets disappears from the Image Selector icon and the icon is again constructed from the image frame that was acquired at 2.5 s after radiography has started.

Keep NumLock set (the default) on the bedside console.


### 4.6.5 Filtering a Reference Image

To filter the Image Selector so only icons for images that were added as reference images, follow this procedure.

- 1 Press the IVR NEO [Reference Image] button or the SMART Touch [Filter Ref Image] button. (LED of IVR NEO turns on)
- 2 Alternatively, click the [Image Filtering] button on the bottom GUI bar and then in the Image Filtering list, choose [Having Reference Image].




The Image Selector is filtered so only the icons of images that were added as reference image.



- 3 To revert to an unfiltered Image Selector, press the IVR NEO [Reference Image] button or the SMART Touch [Filter Ref Image] button again.

IVR NEO of LED turns off. Or click the [Image Filtering] button on the bottom GUI bar and then in the Image Filtering list, choose [All Images].

The Image Selector is no-longer filtered; all icons are shown like this.





















## 4.7 Acquisition Keys and GUI Buttons





The following tables summarize the keyboard keys and GUI buttons available to interact with a study that is open for acquisition. All keyboard keys referenced are in the keyboard edit group between the main keys and the keypad.

A: Acquisition monitor, R: Reference monitor

### Menu Interaction (If Digital Program Menu is not displayed)






| No. | Action   | Keyboard  | GUI   | A | R                     |
|-----|--|---|---|---|-----------------------|
| 1   | Display Fluoro program menu.   |   |    |   | <input type="radio"/> |
| 2   | Display Rad program menu.  |   |    |   | <input type="radio"/> |
| 3   | Display Image Comment menu   |   |    |   | <input type="radio"/> |
| 4   | Play/show next loop/image.   |    |    |   | <input type="radio"/> |
| 5   | Play/show previous loop/image.   |   |   |   | <input type="radio"/> |
| 6   | If playing, pause.   |  or  |    |   | <input type="radio"/> |
| 7   | If paused, show next/previous frame.   |  /   |  /  |   | <input type="radio"/> |
| 8   | If paused, start playing. (Keyboard only: Hold key for more than 2 seconds to begin forward/backward and play) |  /   |  /  |   | <input type="radio"/> |

\*From No.4 to 8 keyboard keys can be used in combination with the [Ctrl] key to perform the same functions on the Acquisition monitor.





| No. | Action   | Keyboard | GUI   | A | R                     |
|-----|--|----------|---|---|-----------------------|
| 1   | Playback speed: increase (right-click)/decrease (left-click) |          |  |   | <input type="radio"/> |
| 2   | Brightness: increase (right-click) / decrease (left-click).  |          |  |   | <input type="radio"/> |
| 3   | Contrast: increase (right-click) / decrease (left-click).    |          |  |   | <input type="radio"/> |
| 4   | Zoom: increase (right-click) / decrease (left-click).        |          |  |   | <input type="radio"/> |

\*These keyboard keys can be used in combination with the [Ctrl] key to perform the same functions on the Acquisition monitor.

### Menu Interaction (If Rad, Fluoro, Comment Menu are displayed)







| No. | Action  | Keyboard   | GUI | A | R                     |
|-----|---|--|-----|---|-----------------------|
| 1   | Highlight next menu option.   |  |     |   | <input type="radio"/> |
| 2   | Highlight previous menu option.   |  |     |   | <input type="radio"/> |
| 3   | Highlight the top menu option.  |  |     |   | <input type="radio"/> |
| 4   | Highlight the bottom menu option.   |  |     |   | <input type="radio"/> |
| 5   | For Comments, select highlighted comment and close menu.<br>For Rad / Fluoro first-level menu, display second-level menu for highlighted item. For Rad / Fluoro second-level menu, select highlighted program and close menu. |  |     |   | <input type="radio"/> |

### Reference Images

| No. | Action   | Keyboard   | GUI | A | R                     |
|-----|--|--|-----|---|-----------------------|
| 1   | Insert Reference image.  |   |     |   | <input type="radio"/> |
| 2   | Delete Reference image.  |  |     |   | <input type="radio"/> |
| 3   | Go to next loop with image selected as Reference and pause on the image.     |  |     |   | <input type="radio"/> |
| 4   | Go to previous loop with image selected as Reference and pause on the image. |  |     |   | <input type="radio"/> |


\*No. 3 and 4 keys can be used in combination with the [Ctrl] key to perform the same functions on the Acquisition monitor.

### 16 Up Select

| No. | Action  | Keyboard   | GUI | A | R                     |
|-----|---|--|-----|---|-----------------------|
| 1   | Move selector box down one row. (When at bottom row, wrap around to top.)           |  |     |   | <input type="radio"/> |
| 2   | Move selector box up on row. (When at top, wrap around to bottom.)                  |  |     |   | <input type="radio"/> |
| 3   | Move selector box to right one column. (When at right column, wrap around to left.) |  |     |   | <input type="radio"/> |
| 4   | Move selector box to left one column. (When at left column, wrap around to right.)  |  |     |   | <input type="radio"/> |
| 5   | Move selector box down four rows. (When at down, wrap around to top.)               |  |     |   | <input type="radio"/> |
| 6   | Move selector box up four rows. (When at top, wrap around to bottom.)               |  |     |   | <input type="radio"/> |

## 4.8 IVR NEO (Either-or SMART Touch)

The default IVR NEO features 16 pre-configured buttons (each with an indicator LED), and a joystick with button. Default button settings are as follows.



4

See "Physicians Tab" P.17-44 for changing the default settings of these buttons.



































### NOTE

The IVR NEO and IVR Shuttle operation is replicated to one another, as they work with each other. For example, when switching the IVR NEO focus from the Reference monitor to the Acquisition monitor with the IVR NEO [Monitor Select] button, the IVR Shuttle operation is activated on the Acquisition monitor as well. Therefore, do not use the IVR NEO and IVR Shuttle at the same time.



















### 4.8.1 IVR NEO Buttons

IVR NEO button and LED indicator functionality is summarized in the following table (see button numbers indicated in above IVR NEO illustration). However, functionality availability varies depending on the type of studies you have chosen, such as cardiac or head, so some buttons may not be activated.




















| Category    | Button Name           | LED State   | Icon  | A                     | R                     |
|-------------|-----------------------|---|---|-----------------------|-----------------------|
| Fluoroscopy | Select Fluoro Program | ON= Fluoro menu shown. Goes OFF once selection made.  |    |                       | <input type="radio"/> |
|             | DUP-Fluoro            | Associated predefined Fluoro mode selected upon selection.  |   | <input type="radio"/> | <input type="radio"/> |
|             | Save Fluoro           | ON= Save available. ([Fluoro]-[Save] buttons on the side menu is also available.)   |    | <input type="radio"/> | <input type="radio"/> |
|             | Save Last Fluoro      | ON= Fluoroscopy images can be saved mode.   |    | <input type="radio"/> | <input type="radio"/> |
|             | View Fluoro           | Last acquired Fluoro loop shown upon selection.<br>ON= Last acquired Fluoro loop shown.   |    |                       | <input type="radio"/> |
|             | Direct FLUORO Record  | Direct Record mode enabled upon selection.<br>ON= Direct Record Mode.<br>LED goes OFF upon completion of next Fluoro acquisition. |  |                       | <input type="radio"/> |
| Radiography | Select Rad Program    | ON= Rad menu shown. Goes OFF once selection made.   |  |                       | <input type="radio"/> |
|             | DUP-Rad               | Associated predefined Rad mode selected upon selection.   |   | <input type="radio"/> | <input type="radio"/> |
|             | Change Planes         | (Biplane only)<br>Change the plane in order of Frontal, Bi and Lateral.   |  | <input type="radio"/> | <input type="radio"/> |
|             | Bi-plane              | (Biplane only)<br>Change the plane to Bi-plane.   |  | <input type="radio"/> | <input type="radio"/> |
|             | Frontal               | (Biplane only)<br>Change the plane to Frontal.  |  | <input type="radio"/> | <input type="radio"/> |
|             | Lateral               | (Biplane only)<br>Change the plane to Lateral.  |  | <input type="radio"/> | <input type="radio"/> |
| MAP         | DSA-MAP               | ON= DSA-MAP mode.   |  | <input type="radio"/> | <input type="radio"/> |
|             | DSA-MAP (Live)        | ON= DSA-MAP Live mode.  |  | <input type="radio"/> | <input type="radio"/> |
|             | DSA-MAP (Sub)         | ON= DSA-MAP Sub mode.   |  | <input type="radio"/> | <input type="radio"/> |

| Category | Button Name            | LED State   | Icon  | A                     | R                     |
|----------|------------------------|---|---|-----------------------|-----------------------|
| MAP      | FluoroMAP              | ON= FluoroMAP Mask Creation mode.   |    | <input type="radio"/> | <input type="radio"/> |
|          | FluoroMAP (Live)       | ON= FluoroMAP Live mode.  |    | <input type="radio"/> | <input type="radio"/> |
|          | FluoroMAP (Sub)        | ON= FluoroMAP Sub mode.   |    | <input type="radio"/> | <input type="radio"/> |
|          | BlankMAP               | Register the Blank mask.  |    | <input type="radio"/> | <input type="radio"/> |
|          | MAP mode               | ON= MAP mode.   |    | <input type="radio"/> | <input type="radio"/> |
|          | LIVE mode              | ON= Live mode.  |    | <input type="radio"/> | <input type="radio"/> |
|          | Reset BG Mask          | Reset the background mask.  |    | <input type="radio"/> | <input type="radio"/> |
|          | TRANS                  | ON= IVR NEO joystick changes MAP Power.   |    | <input type="radio"/> | <input type="radio"/> |
|          | Sketch Display         | ON= Sketch Display mode.  |   | <input type="radio"/> | <input type="radio"/> |
|          | Sketch Edition         | ON= Sketch Edition window is displayed, and the Joystick becomes available.       |  | <input type="radio"/> | <input type="radio"/> |
|          | Contour                | Display the contour enhanced image.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Mask Region            | ROI of SIMAP is displayed on the image and the joystick is in the operation mode. |  | <input type="radio"/> | <input type="radio"/> |
| DSA      | LIVE                   | ON= Display DSA Live image.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Peak Hold              | Display the peak hold image.  |  | <input type="radio"/> | <input type="radio"/> |
|          | Peak Hold (All Frames) | Display the peak hold image which used all frames.                                |  | <input type="radio"/> | <input type="radio"/> |
|          | Flex APS (Option)      | ON= Display the position offset image.  |  | <input type="radio"/> | <input type="radio"/> |
| Actions  | X-ray OFF              | ON= Ban the X-ray exposure.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Start/Stop/Close Timer | ON= Timer shown on the Acquisition monitor.                                       |  | <input type="radio"/> | <input type="radio"/> |
|          | C-arm Reposition       | Send the angle of displayed image to the C-arm.                                   |  | <input type="radio"/> | <input type="radio"/> |



## 4 Image Acquisition

| Category | Button Name             | LED State   | Icon  | A                     | R                     |
|----------|-------------------------|---|---|-----------------------|-----------------------|
| Actions  | Add Reference Image     | ON= Add Reference Image available.  |    | <input type="radio"/> | <input type="radio"/> |
|          | Delete reference Image  | Reference image deleted upon selection.<br>ON= Delete Reference Image available.                |    | <input type="radio"/> | <input type="radio"/> |
|          | Save Still Image        | ON= Save available.<br>Blank once upon selecting.   |    | <input type="radio"/> | <input type="radio"/> |
|          | Add Print               | Image is added for printing upon selection.<br>ON= Add Print available.                         |    | <input type="radio"/> | <input type="radio"/> |
|          | Monitor Select          | ON= Acquisition monitor selected.<br>OFF= Reference.  |    | <input type="radio"/> | <input type="radio"/> |
|          | Lateral Operation       | (Bi-plane only)<br>Enable to operate from ACQ-L when "Edit Sketch" and "MAP Power" are changed. |    | <input type="radio"/> | <input type="radio"/> |
|          | Reference Mode          | ON= Reference Mode.   |    | <input type="radio"/> | <input type="radio"/> |
|          | Ref Sync                | ON= RefSync mode.   |   | <input type="radio"/> | <input type="radio"/> |
|          | PCI View Mode           | ON= PCI View mode.  |  | <input type="radio"/> | <input type="radio"/> |
|          | Selecting the Last Loop | Select the last loop on the image selector.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Table Reposition        | ON= Send the position of displayed image to the tabletop.                                       |  | <input type="radio"/> | <input type="radio"/> |
|          | Fluoro Alarm            | Cancel fluoro alarm sound.  |  | <input type="radio"/> | <input type="radio"/> |
| Display  | Filter by C-arm         | ON= Images filtered to match C-arm position.  |  | <input type="radio"/> | <input type="radio"/> |
|          | Reference Image         | Search and display the added reference images.  |  | <input type="radio"/> | <input type="radio"/> |
|          | 16-up Selector          | ON= 16-up Select mode.<br>Goes OFF once selection made.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Single/Split Display    | IVR NEO focus= Reference.<br>Reference split.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Left/Right Focus        | ON= Right half of split Reference monitor selected.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Cine Area Zoom          | ON= Zoom the image in the Examination room.   |  | <input type="radio"/> | <input type="radio"/> |



| Category | Button Name               | LED State  | Icon  | A                     | R                     |
|----------|---------------------------|--|---|-----------------------|-----------------------|
| Display  | Zoom In                   | Zoom level increases slightly upon selection.<br>ON= Zoom In available.    |    | <input type="radio"/> | <input type="radio"/> |
|          | Zoom Out                  | Zoom level decreased slightly upon selection.<br>ON= Zoom out available.   |    | <input type="radio"/> | <input type="radio"/> |
|          | H-Flip                    | Flip an acquisition image horizontally.                                    |    | <input type="radio"/> | <input type="radio"/> |
|          | V-Flip                    | Flip an acquisition image vertically.                                      |    | <input type="radio"/> | <input type="radio"/> |
| Play     | Play Forward              | Loop plays forward upon selection.<br>ON= Play Forward available.          |    | <input type="radio"/> | <input type="radio"/> |
|          | Play Backward             | Loop plays backward upon selection.<br>ON= Play Backward available.        |    | <input type="radio"/> | <input type="radio"/> |
|          | Play/Sync Mode            | ON= Play Single Loop with Synchronize mode.<br>OFF= Play Single Loop mode. |    | <input type="radio"/> | <input type="radio"/> |
| Analysis | LV Mode                   | LV menu shown upon selection.<br>ON= LV menu shown.                        |   | <input type="radio"/> | <input type="radio"/> |
|          | QCA Mode                  | QCA menu shown upon selection.<br>ON= QCA menu shown.                      |  | <input type="radio"/> | <input type="radio"/> |
|          | Select Image Comment      | ON= Comment menu shown. Goes OFF once selection made.                      |  | <input type="radio"/> | <input type="radio"/> |
|          | Add annotation's Text     | Add text tool selected upon selection.                                     |  | <input type="radio"/> | <input type="radio"/> |
|          | Add annotation's Line     | Line tool added at mouse cursor location upon selection.                   |  | <input type="radio"/> | <input type="radio"/> |
|          | Add annotation's Curve    | Curve tool added at mouse cursor location upon selection.                  |  | <input type="radio"/> | <input type="radio"/> |
|          | Add annotation's Arrow    | Arrow added at mouse cursor location upon selection.                       |  | <input type="radio"/> | <input type="radio"/> |
|          | Add annotation's Caliper  | Caliper tool added at mouse cursor location upon selection.                |  | <input type="radio"/> | <input type="radio"/> |
|          | Add annotation's Analysis | Analysis results added upon selection.                                     |  | <input type="radio"/> | <input type="radio"/> |
|          | Add annotation's Study    | Study information added upon selection.                                    |  | <input type="radio"/> | <input type="radio"/> |
|          | Clear last annotation     | Last annotation deleted upon selection.                                    |  | <input type="radio"/> | <input type="radio"/> |
|          | Clear all annotations     | All annotations deleted upon selection.                                    |  | <input type="radio"/> | <input type="radio"/> |

## 4 Image Acquisition

| Category | Button Name                | LED State                                    | Icon  | A | R                     |
|----------|----------------------------|--|---|---|-----------------------|
| Analysis | Show Annotation properties | Annotations properties shown upon selection. |  |   | <input type="radio"/> |
|          | Add annotation's Pointer   | Pointer tool selected upon selection.        |  |   | <input type="radio"/> |

A: relevant to Acquisition monitor, R: relevant to Reference monitor



Change DUP or radiography mode during fluoroscopy by shortcut keys and IVR NEO buttons, after completing current fluoroscopy.



Some shortcut keys are available at the same time with [AcqSync] button on the REF-PC side menu.

## 4.8.2 IVR NEO Joystick

The IVR NEO joystick can be used to playback images and perform other actions as follows.

| No | Category  | Initial GUI State   | Joystick Function   | A | R |
|----|-----------|---|---|---|---|
| 1  | Menu      | Fluoro Program menu displayed   | Up = move to previous menu item.<br>Down = move to next menu item.  |   | ○ |
|    |           | Rad Program menu displayed  | Left = move to first menu item.<br>Right = move to last menu item.  |   | ○ |
|    |           | Image Comment menu displayed  | Click joystick button = select highlighted menu item.   |   | ○ |
| 2  | Playback  | Cyclic Mode:<br>(Playback in forward direction at acquisition speed.) | Up / Down = Play previous / next loop. In reference mode, display the previous / next reference image.<br>Click joystick button = Switch to Still mode, pausing on image shown when clicked button.   | ○ | ○ |
|    |           | Still Mode:<br>(Paused on still image in loop.)                       | Up / Down = Display previous / next loop. If a loop contains a reference image, display the reference image. In reference mode, display the previous / next reference image only.<br><br>Left / Right then release (below threshold) = Display previous frame (to left) or next frame (to right) and then pause.<br><br>Left / Right and hold (above threshold) = Play frame-by-frame backward (to left) / forward (to right) at speed (10% to 200%) set by joystick volume. Release joystick to Pause.<br><br>Click joystick button = Switch to Cyclic mode and play loop forward at acquired speed. | ○ | ○ |
| 3  | Selection | 16-up Selector window displayed                                       | Left, Right, Up, Down = Move selection box around screen.   |   | ○ |
|    |           |   | Up / Down at screen top/bottom = Scroll Selection screen if more than 16 images.  |   | ○ |
|    |           |   | Click joystick button = Switch to Cyclic Mode and play selected loop (in single-image Image Viewer window) at acquired speed.   |   | ○ |
| 4  | ROI       | ROI Move Mode   | Left, Right, Up, Down = Move ROI in parallel.   |   | ○ |
|    |           |   | Click joystick button = Switch to ROI Size Modification Mode.   |   | ○ |
|    |           | ROI Size Modification Mode  | Left, Right, Up, Down = Modify ROI size by operating the peak of bottom-right of ROI as fixing the peak of upper-left of ROI.   |   | ○ |
| 5  | MAP Power | MAP Power Modification Mode   | Left= Decrease MAP power.<br>Right= Increase MAP power.   | ○ |   |

A: relevant to Acquisition monitor, R: relevant to Reference monitor. Playback of the image is common in both normal and split display.

## 4.9 Keyboard ShortCuts

Keyboard function keys have equivalent operation as IVR NEO buttons. ([Monitor Select] and [Reference Mode] functions are supported only for IVR NEO)

Move the mouse to Acquisition monitor to operate on Acquisition monitor.


 "IVR NEO Buttons" P.4-66.

The following functions are set to IVR NEO and Shortcut Key as default. It is enable to set functions by holding down the left Shift key and press a function key (F1~F12).

| IVR NEO |                        | Shortcut Key                     |                        |
|---------|------------------------|----------------------------------|------------------------|
| Key     | Description            | Key                              | Description            |
| B1      | Select Fluoro Program  | F1                               | Select Fluoro Program  |
| B2      | Monitor Select         | F2                               | Select Rad Program     |
| B3      | Select Rad Program     | F3                               | Live                   |
| B4      | Live                   | F4                               | Filter by C-arm        |
| B5      | Filter by C-arm        | F5                               | C-arm Reposition       |
| B6      | C-arm Reposition       | F6                               | Start/Stop/Close Timer |
| B7      | Start/Stop/Close Timer | F7                               | Cine Area Zoom         |
| B8      | Cine Area Zoom         | F8                               | Select Image Comment   |
| B9      | Select Image Comment   | F9                               | 16-up Selector         |
| B10     | 16-up Selector         | F10                              | Save Fluoro            |
| B11     | Save Fluoro            | F11                              | Save Still Image       |
| B12     | Save Still Image       | F12                              | Add Reference Image    |
| B13     | Add Reference Image    | Left Shift+F1~<br>Left Shift+F12 | Not configured         |
| B14     | Reference Mode         |                                  |                        |
| B15     | Single/Split Display   |                                  |                        |
| B16     | Left/Right Focus       |                                  |                        |

## 4.10 IVR Shuttle (Option) (Either-or SMART Touch)
















The IVR Shuttle features 17 buttons (each with an indicator LED) and 14 controls, including a dual-control shuttle. Button and control settings are as follows.




The IVR NEO and IVR Shuttle operation is replicated to one another, as they work with each other. For example, when switching the IVR NEO focus from the Reference monitor to the Acquisition monitor with the IVR NEO [Monitor Select] button, the IVR Shuttle operation is activated on the Acquisition monitor as well. Therefore, do not use the IVR NEO and IVR Shuttle at the same time.

## 4.10.1 Functions of IVR Shuttle Buttons and Controls

## 1 Buttons

| Category    | Button Name             | Description  | Icon   | A | R |
|-------------|-------------------------|--|--|---|---|
| Acquisition | Select Fluoro Program   | ON= Fluoro menu is displayed. Scroll up and down the menu selector (2) to move to the menu items. Press the menu selector to select the highlighted menu.                |  Fluo     |   | ○ |
|             | Select Rad Program      | ON= Rad menu is displayed. Scroll up and down the menu selector (2) to move to the menu items. Press the menu selector to select the highlighted menu.                   |  Rad      |   | ○ |
|             | Add Reference Image     | ON= Add Reference Image available.   |           | ○ | ○ |
|             | Save Fluoro             | ON= Save available. Fluoroscopy loop can be saved in buffer.   |           | ○ | ○ |
|             | Start/Stop/ Close Timer | Press once to display the timer on acquisition monitor to start. Press again to stop and press third time to hide the timer.   |           | ○ |   |
| Monitor     | Monitor Select          | ON= Acquisition monitor is selected.<br>OFF= Reference monitor is selected.  |         | ○ | ○ |
|             | Zoom                    | ON= Zoom is available.<br>Scroll up and down the menu selector(2) to change the magnification.   |  ZOOM   | ○ | ○ |
|             | Live                    | ON= DSA image is displayed.<br>OFF= LIVE image is displayed.   | LIVE   | ○ | ○ |
|             | REMASK                  | ON= Create the mask again.   |  REMASK |   | ○ |
|             | Pause                   | ON= LIVE image is paused.  |         | ○ | ○ |
|             | Select Reference Image  | Select reference image. Display the previous image.  |         | ○ | ○ |
|             |                         | Display the next image.  |         | ○ | ○ |
|             | Select Loop             | ON= Acquisition monitor is selected.<br>OFF= Reference monitor is selected.  |         | ○ | ○ |
|             |                         | Review next loop.  |         | ○ | ○ |
|             | Frontal                 | For Bi-plane system, select a target processing for F or L individually. Adjusting Brightness/ Contrast, modifying Edge Enhancement and Save Still Image are processing. |  F      |   | ○ |
|             | Lateral                 |  |  L      |   | ○ |

| Category | Button Name      | Description         | Icon  | A | R                     |
|----------|------------------|---------------------|---|---|-----------------------|
| Others   | Save Still Image | ON= Save available. |  |   | <input type="radio"/> |

A: relevant to Acquisition monitor, R: relevant to Reference monitor

## 2 Menu



Scroll up and down to move to the menu items. Press the menu selector to select the highlighted menu.


## 3 Change Frame

Turn the internal dial to the right to display the next frame or to the left to display the previous frame.

## 4 Loop Speed Adjust

Turn the external dial to the right to play the loop forward and to the left to play the loop backward. Turn more in either direction to increase speed.

## 5 Brightness/Contrast Rest

Turn the [Brightness] dial to the right to increase image brightness or to the left to decrease brightness (max. 100, min. -100). And turn the [Contrast] dial to the right to increase image contrast or to the left to decrease contrast (max. 2000, min. -100). Press [Reset]  to return to the default brightness and contrast values.

## 6 Edge

Adjust the [Convolution] filter. Turn the [Edge] dial to the right to increase the sharpness of the image and to the left to decrease sharpness (max. 100, min. -100). Press the dial to return to the default value.

## 7 Change Field of View



Press the right button to increase the size of the FOD field of view and left reduce.

## 8 Collimator C-leaf and Peripheral Filter



Use the upper button to open the C-leaf collimator and the lower button to close it.

## 9 Collimator V/H-leaf



Open and close the collimator V/H-leaf. Use the up (open) /down (close) button for the collimator V-leaf and the right (open) / left (close) button for the collimator H-leaf. Press the center button to fully open the H, V and C leaf.

## 10 Compensation Filter

Open, close and rotate the right and left compensation filter.


### 11 Exam Room Switch Monitor

If there are option computers such as 3D workstation, connect display in the examination room.






## 4.11 SMART Touch (Either-or IVR NEO/IVR Shuttle)

SMART Touch is consist of Tab Area, Panel Area and Side Menu. And a joystick is connected to SMART Touch, and menu selection and image control are available.



The following table shows description of each area.

| No. | Function   | Outline   |
|-----|------------|---|
| 1   | Tab Area   | <p>Consist of message area and tabs to switch the panel display. Type of tabs are as follows:<br/>Custom Tab, Program Tab, Tool Tab and Monitor Tab (when connected to SMART Display)</p> <p>Following messages are displayed in the Message Area.<br/>When connected to DAR-9500f system:</p>  <p>Before start new study:</p>  <p>During study:</p>  |
| 2   | Panel Area | Contents of display will change depending on the selected tab.  |

## 4 Image Acquisition

---

| No. | Function  | Outline   |
|-----|-----------|---|
| 3   | Side Menu | Consist of [X-ray OFF], [Option] buttons, and 5 of function buttons. And if fluoroscopy time exceeds the setting fluoroscopy time, [Fluoro Alarm] button will be displayed. |



Enable to connect up to 3 SMART Touch.

As each SMART Touch synchronized, operation of SMART Touch reflects to other SMART Touch. (For example, if change the target monitor from REF monitor to ACQ monitor with [Monitor Select] button of SMART Touch which s placed in the control room, SMART Touch which is placed in the examination room will also control the ACQ monitor.) So, do not use several SMART Touch at the same time.




Function buttons are configurable individually for each SMART Touch.

### 4.11.1 Custom Panel

Custom panel is displayed shortly after the system startup and touching the Custom Tab. An optional function can be registered for each physician on the Custom panel.

The following figure shows default settings.



Refer to "4.11.6 SMART Touch Buttons" P.4-100 for each function.

#### Option Button



Unable to use [Option] button during study.

Modify the following settings from the Option window.


- Assignment of buttons on the Custom Panel and Side Menu
- Layout of Custom Panel
- Color Theme
- Volume

Configuration is available for each physician on the Option window.



And if install several SMART Touch, enable to set for each SMART Touch appropriate to usage and operator.




## 4 Image Acquisition


Tab Area is consist of [General option] and [Custom tab option].





| No.                      | Area           | Outline                          |
|--------------------------|----------------|----------------------------------|
| <b>Custom Tab Option</b> |                |                                  |
| 1                        | Physician Name | Select a physician to customize. |

| No. | Area   | Outline  |
|-----|--------|--|
| 2   | Target | <p>Select a target to customize from Panel Area and Sidebar.</p> <p>Panel Area:</p>  <p>Physician Name: Default<br/>Target: Panel (selected), Sidebar<br/>Layout: 2x2, 4x3, 5x4 (selected)</p> <p>Layout preview</p> <p>Actions:</p> <ul style="list-style-type: none"><li>X-Ray OFF</li><li>Stopwatch</li><li>C-Arm Reposition</li><li>Table Reposition</li><li>Dose Reset</li></ul> <p>1/4</p> <p>Sidebar:</p>  <p>Physician Name: Default<br/>Target: Panel, Sidebar (selected)<br/>Layout: 2x2, 4x3, 5x4</p> <p>Layout preview</p> <p>Fluoro</p> <ul style="list-style-type: none"><li>Select Fluoro Program</li><li>Save Fluoro</li><li>Save Fluoro Last 1F</li><li>View Fluoro</li><li>Direct Fluoro Rec</li></ul> <p>1/1</p> |

| No. | Area           | Outline   |
|-----|----------------|---|
| 3   | Layout         | <p>Select the number of buttons to display on the Custom Panel from [2x2], [4x3] and [5x4]<br/>2x2:..</p>  <p>4x3:</p>  <p>5x4:</p>  |
| 4   | Layout preview | Display the buttons assigned to the Custom Panel and Sidebar.   |
| 5   | Function       | Display functions that are to be assigned.  |




| No.            | Area        | Outline  |
|----------------|-------------|--|
| General Option |             |  |
| 1              | Color Theme | <p>Modify the display color.</p> <p>Dark:</p>  |

| No. | Area   | Outline  |
|-----|--------|--|
|     |        | <p>Light:</p>  |
| 2   | Volume | Modify the volume of Touch sound.  |

### ■ Modify Assignment of Buttons


Assign optional buttons to Custom Panel and Side Menu for each physician and SMART Touch.

#### 1 Touch [Option].







- 2 Select customized physician from [Physician Name].




- 3 Select [Panel] or [Sidebar] from [Target]. Display of Layout preview is changed.




### 4 Select a category from [Function].




### 5 Select a button to be modified from [Layout preview].



6 And then select a button to be assigned.




The button is modified.





### Change of Button Colors


Change the color of buttons on the Custom Panel. If display color is Dark, enable to set the following colors.




If display color is Light, enable to set the following colors.




- 1 Select a button to be modified from [Layout preview].



2 Select the color of button.




The color of button will be changed.



### 4.11.2 Program Panel

Program Panel is displayed when touching the Program tab. There are buttons to switch Digital User Program (DUP), Radiography Plane, and FOV.




#### Switch RAD DUP

Touch RAD menu to display a pull-down menu and touch desired DUP.

#### Switch RAD Plane (In case of Bi-plane system)

Touch the plane to perform radiography.




#### Switch Fluoro DUP

Touch Fluoro menu to display a pull-down menu and touch desired DUP.


#### Switch FOV

Touch FOV to be used.



### 4.11.3 Tool Panel

Tool panel will display when touch the Tool Tab. There are pre-assigned buttons on the tool tab.



The following table shows the registered functions.

| Function                                   | Outline  |
|--|--|
| <b>Tool Area</b>                           |  |
| Stopwatch                                  | Timer shown on the acquisition monitor.        |
| Photo Capture                              | Save still image.                              |
| Save Fluoro                                | Save fluoroscopy image.                        |
| <b>Reposition Area</b>                     |  |
| C-Arm Reposition                           | Move the C-arm to position of the image.       |
| Table Reposition<br>(Combined with KS-100) | Move the table to position of the image.       |
| <b>Reference Image Area</b>                |  |
| Add Ref Image                              | Register reference image.                      |
| Delete Ref Image                           | Delete reference image.                        |
| <b>Image Selector Area</b>                 |  |
| Filter by C-Arm                            | Filter the image based on the C-arm.           |
| Filter Ref Image                           | Filter the image based on the reference image. |
| Select Last Loop                           | Select the last acquisition image.             |
| Play Sync/Play                             | Switch [Sync] and [Play].                      |

| Function                    | Outline   |
|-----------------------------|---|
| <b>Monitor Control Area</b> |   |
| ACQ Sync                    | (For Bi-plane system only)<br>Synchronize the operation of ACQ-F monitor and ACQ-L monitor. |
| REF Sync                    | (For Bi-plane system only)<br>Synchronize the operation of REF-F monitor and REF-L monitor. |
| ACQ                         | Operation target is ACQ monitor.  |
| REF                         | Operation target is REF monitor.  |
| Frontal                     | Operation target is Frontal.  |
| Lateral                     | (For Bi-plane system only)<br>Operation target is Lateral.                                  |




#### 4.11.4 Monitor Panel (Only when SMART Display is connected)



Be sure to include acquisition monitor in the preset.

Monitor panel will display when touch the Monitor Tab. Select SMART Display preset and switch video input.




The following table shows description of each area.


| No. | Area           | Outline  |
|-----|----------------|--|
| 1   | Preset List    | Select preset of SMART Display. Selected preset is displayed on the preset preview.      |
| 2   | Preset Preview | Display preset of SMART Display. Display video input on each area in the preset preview. |
| 3   | Input List     | List video input that can be displayed on the SMART Display.                             |

### Change of Preset

- 1 Select preset from the preset list.  
Touch an arrow key to change the page of preset list.




- 2 Select preset to display selected preset on the preset preview.




### Change of Segment

- 1 Select segment to change the display position on the preset preview.



2 Select segment to be displayed.




Segment is changed.




4

### Change of Video Input


- 1 Select video input to be changed on the preset preview.



- 2 Select video input to be displayed from the input list.




Video input is changed.




### 4.11.5 Side Menu

There are [X-ray OFF] button and function buttons on the side menu. And if exceed setting fluoroscopy time, [Fluoro Alarm] stop button display.




#### ■ X-ray OFF

Touch [X-ray OFF] button to stop X-ray radiography and fluoroscopy. X-ray OFF icon turns to red if it is selected. Touch the button again to enable X-ray radiography and fluoroscopy.



## Function Button

The following figure shows default setting.



Refer to  ["4.11.6 SMART Touch Buttons" P.4-100](#) for each function.

## Stop Fluoro Alarm














If exceed setting fluoroscopy time, [Fluoro Alarm] stop button will display. Touch [Fluoro Alarm] stop button to stop fluoroscopy alarm.






















### 4.11.6 SMART Touch Buttons

SMART Touch buttons are summarized in the following table.




















Color of button turns green when it is available. And button turns grayout when it is not available.

















| Category    | Button Name           | LED State  | Icon  | A                     | R                     |
|-------------|-----------------------|--|---|-----------------------|-----------------------|
| Fluoroscopy | Select Fluoro Program | Display fluoro menu.   |    |                       | <input type="radio"/> |
|             | DUP-Fluoro            | Associated predefined Fluoro mode selected upon selection.   |   | <input type="radio"/> | <input type="radio"/> |
|             | Save Fluoro           | Record the latest fluoro loop. ([Fluoro]-[Save] buttons on the side menu is also available.)                 |    | <input type="radio"/> | <input type="radio"/> |
|             | Save Last Fluoro      | Record the latest frame of fluoro loop.  |    | <input type="radio"/> | <input type="radio"/> |
|             | View Fluoro           | Last acquired Fluoro loop shown upon selection.  |    |                       | <input type="radio"/> |
|             | Direct Fluoro Rec     | Save at the same time with fluoroscopy acquisition.<br>It becomes unavailable after fluoroscopy acquisition. |   |                       | <input type="radio"/> |
| Radiography | Select Rad Program    | ON= Rad menu shown. Goes OFF once selection made.  |  |                       | <input type="radio"/> |
|             | DUP-Rad               | Associated predefined Rad mode selected upon selection.  |   | <input type="radio"/> | <input type="radio"/> |
|             | Toggle Rad Mode       | (Biplane only)<br>Change the plane in order of Frontal, Bi and Lateral.                                      |  | <input type="radio"/> | <input type="radio"/> |
|             | Rad Biplane           | (Biplane only)<br>Change the plane to Bi-plane.  |  | <input type="radio"/> | <input type="radio"/> |
|             | Rad Frontal           | (Biplane only)<br>Change the plane to Frontal.   |  | <input type="radio"/> | <input type="radio"/> |
|             | Rad Lateral           | (Biplane only)<br>Change the plane to Lateral.   |  | <input type="radio"/> | <input type="radio"/> |
| MAP         | DSA-MAP               | DSA-MAP mode.  |  | <input type="radio"/> | <input type="radio"/> |
|             | DSA-MAP (Live)        | DSA-MAP Live mode.   |  | <input type="radio"/> | <input type="radio"/> |
|             | DSA-MAP (Sub)         | DSA-MAP Sub mode.  |  | <input type="radio"/> | <input type="radio"/> |



| Category | Button Name            | LED State   | Icon  | A                     | R                     |
|----------|------------------------|---|---|-----------------------|-----------------------|
| MAP      | FluoroMAP              | FluoroMAP Mask Creation mode.   |    | <input type="radio"/> | <input type="radio"/> |
|          | FluoroMAP (Live)       | FluoroMAP Live mode.  |    | <input type="radio"/> | <input type="radio"/> |
|          | FluoroMAP (Sub)        | FluoroMAP Sub mode.   |    | <input type="radio"/> | <input type="radio"/> |
|          | BlankMAP               | Register the Blank mask.  |    | <input type="radio"/> | <input type="radio"/> |
|          | MAP mode               | MAP mode.   |    | <input type="radio"/> | <input type="radio"/> |
|          | LIVE mode              | Live mode.  |    | <input type="radio"/> | <input type="radio"/> |
|          | Reset BG Mask          | Reset the background mask.  |    | <input type="radio"/> | <input type="radio"/> |
|          | Vessel Power           | SMART Touch joystick changes Vessel Power.  |    | <input type="radio"/> | <input type="radio"/> |
|          | Sketch Display         | Sketch Display mode.  |   | <input type="radio"/> | <input type="radio"/> |
|          | Sketch Edit            | Sketch Edition window is displayed, and the Joystick becomes available.           |  | <input type="radio"/> | <input type="radio"/> |
|          | Contour                | Display the contour enhanced image.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Mask Region            | ROI of SIMAP is displayed on the image and the joystick is in the operation mode. |  | <input type="radio"/> | <input type="radio"/> |
|          |                        |   |   |                       |                       |
| DSA      | LIVE                   | Display DSA Live image.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Peak Hold              | Display the peak hold image.  |  | <input type="radio"/> | <input type="radio"/> |
|          | Peak Hold (All Frames) | Display the peak hold image which used all frames.                                |  | <input type="radio"/> | <input type="radio"/> |
|          | Flex-APS (Option)      | Display the offset position image.  |  | <input type="radio"/> | <input type="radio"/> |
| Actions  | X-ray OFF              | Ban the X-ray exposure.   |  | <input type="radio"/> | <input type="radio"/> |
|          | Timer                  | Timer shown on the Acquisition monitor.   |  | <input type="radio"/> | <input type="radio"/> |
|          | C-Arm Reposition       | Send the angle of displayed image to the C-arm.                                   |  | <input type="radio"/> | <input type="radio"/> |

## 4 Image Acquisition

| Category | Button Name            | LED State   | Icon  | A                     | R                     |
|----------|------------------------|---|---|-----------------------|-----------------------|
| Actions  | Table Reposition       | (When combined with KS-100)<br>Send the position of displayed image to the tabletop.            | <br>Table Reposition       |                       | <input type="radio"/> |
|          | Add Reference Image    | Add Reference Image available.  | <br>Add Ref Image          | <input type="radio"/> | <input type="radio"/> |
|          | Delete reference Image | Reference image deleted upon selection.   | <br>Delete Ref Image       | <input type="radio"/> | <input type="radio"/> |
|          | Photo Capture          | Save still image (must be paused).  | <br>Photo Capture          |                       | <input type="radio"/> |
|          | Add Print              | Image is added for printing upon selection.   | <br>Add Print              |                       | <input type="radio"/> |
|          | Monitor Select         | Select the monitor operating with SMART Touch.  | <br>REF<br>REF             | <input type="radio"/> | <input type="radio"/> |
|          | Lateral Operation      | (Bi-plane only)<br>Enable to operate from ACQ-L when "Edit Sketch" and "MAP Power" are changed. | <br>Frontal<br>Frontal     | <input type="radio"/> |                       |
|          | Reference Mode         | Reference Mode.   | <br>Reference Mode        | <input type="radio"/> | <input type="radio"/> |
|          | REF Sync               | (Bi-plane only)<br>RefSync mode.  | <br>REF<br>REF Sync      |                       | <input type="radio"/> |
|          | ACQ Sync               | (Bi-plane only)<br>RefSync mode.  | <br>ACQ<br>ACQ Sync      |                       | <input type="radio"/> |
|          | PCI View               | PCI View mode.  | <br>PCI View<br>PCI View | <input type="radio"/> |                       |
|          | Select Last Loop       | Select the last loop on the image selector.   | <br>Select Last Loop     | <input type="radio"/> | <input type="radio"/> |
|          | Fluoro Alarm           | Cancel fluoro alarm sound.  | <br>Fluoro Alarm         | <input type="radio"/> | <input type="radio"/> |
| Display  | Filter by C-Arm        | Images filtered to match C-arm position.  | <br>Filter by C-Arm      |                       | <input type="radio"/> |
|          | Filter Ref Image       | Search and display the added reference images.  | <br>Filter Ref Image     | <input type="radio"/> | <input type="radio"/> |
|          | 16-up Selector         | Display in 16-up Select mode.   | <br>16 up selector       |                       | <input type="radio"/> |
|          | Split Display          | Split review monitor in vertically.   | <br>Split Display        |                       | <input type="radio"/> |
|          | Left/Right Focus       | Click a half (left or right) of split image on the Reference monitor.                           | <br>L/R Focus            |                       | <input type="radio"/> |
|          | Cine Area Zoom         | Zoom the images only in the Examination room.   | <br>Cine Area Zoom       | <input type="radio"/> | <input type="radio"/> |

| Category   | Button Name          | LED State   | Icon   | A                     | R                     |
|------------|----------------------|---|--|-----------------------|-----------------------|
| Display    | Zoom In              | Zoom level increases slightly upon selection.     |  Zoom In          | <input type="radio"/> | <input type="radio"/> |
|            | Zoom Out             | Zoom level decreased slightly upon selection.     |  Zoom Out         | <input type="radio"/> | <input type="radio"/> |
|            | H-Flip               | Flip an acquisition image horizontally.           |  H-Flip           | <input type="radio"/> | <input type="radio"/> |
|            | V-Flip               | Flip an acquisition image vertically.             |  V-Flip           | <input type="radio"/> | <input type="radio"/> |
| Play       | Play/Pause           | Switch Play/Pause.                                |  Play             | <input type="radio"/> | <input type="radio"/> |
|            |                      |   |  Pause            | <input type="radio"/> | <input type="radio"/> |
|            | Play Sync            | ON= Play Single Loop with Synchronize mode.       |  Play Sync        | <input type="radio"/> | <input type="radio"/> |
|            | Play Forward         | Loop plays forward upon selection.                |  Play Forward     | <input type="radio"/> | <input type="radio"/> |
|            | Play Backward        | Loop plays backward upon selection.               |  Play Backward   | <input type="radio"/> | <input type="radio"/> |
|            | Previous Loop        | Play/show previous loop/image.                    |  Previous Loop  | <input type="radio"/> | <input type="radio"/> |
|            | Previous Frame       | If paused, show previous frame.                   |  Previous Frame | <input type="radio"/> | <input type="radio"/> |
|            | Next Loop            | Play/show next loop/image.                        |  Next Loop      | <input type="radio"/> | <input type="radio"/> |
|            | Next Frame           | If paused, show next frame.                       |  Next Frame     | <input type="radio"/> | <input type="radio"/> |
| Analysis   | LV Mode              | LV menu shown upon selection.                     |  LV Mode        | <input type="radio"/> | <input type="radio"/> |
|            | QCA Mode             | QCA menu shown upon selection.                    |  QCA Mode       | <input type="radio"/> | <input type="radio"/> |
| Annotation | Select Image Comment | Comment menu shown. Goes OFF once selection made. |  Select Comment | <input type="radio"/> | <input type="radio"/> |

A: relevant to Acquisition monitor, R: relevant to Reference monitor

### 4.11.7 SMART Touch Joystick

The SMART Touch joystick can be used to playback images and perform other actions as follows.

| No | Category  | Initial GUI State   | Joystick Function   | A | R |
|----|-----------|---|---|---|---|
| 1  | Menu      | Fluoro Program menu displayed   | Up = move to previous menu item.<br>Down = move to next menu item.  |   | ○ |
|    |           | Rad Program menu displayed  | Left = move to first menu item.<br>Right = move to last menu item.<br>Click joystick button = select highlighted menu item.   |   | ○ |
|    |           | Image Comment menu displayed  |   |   | ○ |
| 2  | Playback  | Cyclic Mode:<br>(Playback in forward direction at acquisition speed.) | Up / Down = Play previous / next loop. In reference mode, display the previous / next reference image.<br>Click joystick button = Switch to Still mode, pausing on image shown when clicked button.   | ○ | ○ |
|    |           | Still Mode:<br>(Paused on still image in loop.)                       | Up / Down = Display previous / next loop. If a loop contains a reference image, display the reference image. In reference mode, display the previous / next reference image only.<br><br>Left / Right then release (below threshold) = Display previous frame (to left) or next frame (to right) and then pause.<br><br>Left / Right and hold (above threshold) = Play frame-by-frame backward (to left) / forward (to right) at speed (10% to 200%) set by joystick volume. Release joystick to Pause.<br><br>Click joystick button = Switch to Cyclic mode and play loop forward at acquired speed. | ○ | ○ |
| 3  | Selection | 16-up Selector window displayed                                       | Left, Right, Up, Down = Move selection box around screen.   |   | ○ |
|    |           |   | Up / Down at screen top/bottom = Scroll Selection screen if more than 16 images.  |   | ○ |
|    |           |   | Click joystick button = Switch to Cyclic Mode and play selected loop (in single-image Image Viewer window) at acquired speed.   |   | ○ |
| 4  | ROI       | ROI Move Mode   | Left, Right, Up, Down = Move ROI in parallel.   |   | ○ |
|    |           |   | Click joystick button = Switch to ROI Size Modification Mode.   |   | ○ |
|    |           | ROI Size Modification Mode  | Left, Right, Up, Down = Modify ROI size by operating the peak of bottom-right of ROI as fixing the peak of upper-left of ROI.   |   | ○ |
| 5  | MAP Power | MAP Power Modification Mode   | Left= Decrease MAP power.<br>Right= Increase MAP power.   | ○ |   |


A: relevant to Acquisition monitor, R: relevant to Reference monitor. Playback of the image is common in both normal and split display.

## 4.12 Closing the Active Study


Once the diagnostic procedure has been completed, you must close the active study.

Follow this procedure.

- 1 When a new diagnostic procedure will begin immediately, typically for a different patient, click [New] on the side menu.



- 2 The active study is closed and a new study is initiated. When you do not have a new diagnostic procedure to begin immediately, click [Close] instead.



- 3 When MPPS support is enabled, if a Performed Procedure Step manager is used, you will be prompted to select the status of the completed (or discontinued) study like this:

 A screenshot of a dialog box titled 'Modality Performed Procedure Step Information'. It contains a 'Status' dropdown menu set to 'COMPLETED', a 'Protocol Name' text field, a 'Discontinuation Reason' text area, and a 'Performed Protocol Codes' text area. There are right-pointing arrow buttons next to the 'Discontinuation Reason' and 'Performed Protocol Codes' fields. An 'OK' button is at the bottom.

For completed studies, choose [COMPLETED] in the Status list, and optionally, adjust the Performed Protocol codes by clicking the arrow button. Once back at the above window, click [OK] to complete the study.


 A screenshot of a dialog box titled 'Performed Protocol Codes'. It contains a list of protocol codes with checkboxes: 'NECK, T-04001, SCHEME TYPE A, 0.13', 'HAND, T-04002, SCHEME TYPE A, 0.13', 'ABDOMEN, T-04003, SCHEME TYPE A, 0.13', 'ANKLE, T-04004, SCHEME TYPE A, 0.13', and 'FOOT, T-04005, SCHEME TYPE A, 0.13'. 'Modify' and 'Cancel' buttons are at the bottom.

## 4 Image Acquisition

---

For discontinued studies, choose [DISCONTINUED] in the status list and open the Discontinuation Reason list by clicking the arrow button.

Choose the reason and then click [Modify] to return to the above window. Click [OK] to close the discontinues study.



If configured, a Grayscale Softcopy Presentation State is automatically created and saved for all loops just acquired in the study.


Depending on what has already occurred in the study, any number of the following actions may also occur at study closure.

- When MPPS (modality performed procedure step) feature is installed, the PPS (performed procedure step) manager is notified by the MPPS of the study closure and its result, either [Completed] or [Discontinued].
- Reference images are saved as DICOM objects and if configured, are sent to the archive server.
- If configured, media writing (CD/DVD) is automatically started. The media tray opens so that you can insert a blank disc. If your hardware supports it, you may be prompted to choose your media type, CD or DVD, for example. Choose the media type if prompted, insert a blank disc, and close the media tray. If configure, media writing begins automatically once valid blank media is sensed. Otherwise, a pop-up window appears, enabling you to select specific study series to write.

 ["Insert Additional Blank CD/DVDs" P.6-19](#)

## 4.13 Dose Report Image

Dose report image is created automatically after study. Since that the exposure information of study is stored as DICOM image, enable to control exposure of study even if the network server does not support the files in DICOM RDSR format.



This page is intentionally left blank.



# Chapter 5

## Review Basics


This chapter takes you on a tour of the Review process, including finding studies and playing and printing their images.


The key steps of the Review process are introduced here as follows:

### Description

|     |  |     |
|-----|--|-----|
| 5.1 | Startup and Study Selection . . . . .    | 5-2 |
| 5.2 | Viewing Images . . . . .                 | 5-3 |
| 5.3 | Displaying Study Information . . . . .   | 5-6 |
| 5.4 | Creating Image Annotation . . . . .      | 5-6 |
| 5.5 | Saving an Image as a New Object. . . . . | 5-7 |
| 5.6 | Printing Images . . . . .                | 5-7 |
| 5.7 | Using Advanced Features . . . . .        | 5-8 |
| 5.8 | Completing Your Review . . . . .         | 5-8 |


## 5.1 Startup and Study Selection

- 1 In the upper-left corner of the Image Viewer window, click [Open] to display the [Studies Management] window.  
 ["6.1 The Studies Management Window" P.6-2](#)
- 2 Click [Search] to search for all studies performed within the default date range.  
 All matching studies are listed with one row of information per study in descending Date order.



If security is enable, th schedule list is further restricted to only those studies for which the current user is listed as the performing physician


- 3 Click anywhere in the row for the study that you wish to review and then click [View].  
 Alternatively, just double-click any row to view a study.



The first loop of the study begins automatic playback.


## 5.2 Viewing Images

To review a loop other than the first one, click one of the Image-Selector icons. The loop corresponding to the icon begins playing.




5

If there are many loops in the study, Image Selector may continue off the bottom of the window. Right-click (without releasing mouse button) on an icon near the bottom of the window (for example, loop LL7A).




## 5 Review Basics

Drag the mouse upward, scrolling icons for earlier off the top of the window to reveal previously hidden icons.



When you see the icon of a loop that interests you, release the mouse button and then click the icon to begin playing that loop. In this example, we are playing loop LL10A.



5.2.1 Controlling Image Playback

Control image playback with the mouse wheel or via the on-screen Cine Control bar on the bottom border of the Image Viewer window.




For details, including information on keyboard shortcuts and the IVR NEO/IVR Shuttle, see ["7.1 Controlling Image Playback" P.7-2](#).

5.2.2 Adjusting Image Appearance

Adjust image appearance via on-screen [Brightness], [Contrast], [Zoom] and [Filters] button on the side menu, and also [Auto W] button from the [Filters] fly-out menu.

["7.2 Adjusting Image Appearance" P.7-4](#)




5.2.3 Identifying Loops and Still Images

The Image Selector shows small images for both loops and still images. Small icons are super-imposed over the upper-right corner of each small image to indicate whether the small image represents loop or still image.

| No. | Icon | Meaning  |
|-----|------|--|
| 1   |      | A still image created as either an annotated image or reference image. Still images appear in the Image Selector after the loops. The loop number, and in square brackets, the image frame number, are displayed in the upper-left corner of still Image Selector icons. |
| 2   |      | A multi-image rad loop that can be played back as video. The loop number is displayed in the upper-left corner of the Image Selector icon.   |

## 5.3 Displaying Study Information

To display study information including patient demographics, click [Study] on the side menu. The [Study Information] dialog box for the study being reviewed appears.




The image shows the 'Study' menu on the left with the 'Study' icon circled in red. An orange arrow points to the 'Study Information' dialog box on the right. The dialog box contains the following fields:

- Patient**
  - Patient ID: FDDAR9500REV-022
  - Date of Birth: 07/07/1946, Age: 74
  - Patient Name: Prefix (empty), First Name: John, Middle Name (empty), Last Name: Shimadzu, Suffix (empty)
  - Height: inches, Weight: lbs, Sex: Male, Select Orientation
- Study**
  - Study ID: 500001, Accession Number: A000001
  - Modality: RA SR, Referring Physician: Referring First
  - Study Description: (empty)
  - Performing Physician's List: [x] Performing First
  - Operators' List: [x] Operator Second
  - Primary Physician: Performing First
- Buttons: Anonymize, Modify, Close

 "7.4 Displaying Study Information" P.7-10

## 5.4 Creating Image Annotation

To add annotations to a loop image, pause on the image of interest and click [Funcnt]-[ABC Text] on the side menu. You can now add arrows and text and set properties such as color and size for each item you add.





 "7.7 Working with Image Annotations" P.7-17

# 5.5 Saving an Image as a New Object

To save an image, including any annotations, as a new still-image object, click [Funct]-[Save] on the side menu, and then [Yes] to confirm the save.

 "7.8 Saving Images" P.7-20


Alternatively, press  [Save Still Image] button on the IVR NEO/IVR Shuttle/  [Photo Capture] button on the SMART Touch.

# 5.6 Printing Images

To print an image, including any annotations, display the image that you wish to print and then click [Funct] on the side menu. Then click [Print].



The print window appears with an image preview.



Choose the desired printer in the Printer list and click [Print] to print the page as previewed.

 "7.9 Print Images" P.7-21

## 5.7 Using Advanced Features

To use advanced features, consult their respective chapters as follows:

 ["8 Performing Angiographic Analysis"](#)

 ["9 DSA Image Adjustment Tools"](#)

## 5.8 Completing Your Review

When you have finished reviewing a study, click [Close] on the side menu. You can then click [Open] to open the [Studies Management] window and select other studies for review.

This completes the tour of the review process.



# *Chapter 6*

---

## Studies Management


This chapter walks you through the Studies Management process, beginning with a look at its user interface.

### Description

|     |   |     |
|-----|---|-----|
| 6.1 | The Studies Management Window . . . . . | 6-2 |
| 6.2 | Managing Studies . . . . .              | 6-9 |








## 6.1 The Studies Management Window

The [Studies Management] window is displayed by clicking [Open] on the Image Viewer window side menu. It is made up of a columnar Study List, Search Criteria boxes and an Action Bar along the top. Study series and messages are displayed at the bottom.






The above sample screen is made from a composite of images from both monitors. Not all elements appear on both monitors or at the same time. Lines 2, 3, and 4 of the Search Criteria only appear when MPPS support is enabled. ["DICOM Tab" P.17-24](#)

The key components of the Studies Management interface are described here as follows:

| No. | Section                        | Page  |
|-----|--------------------------------|---|
| 1   | "6.1.1 Study Location"         |  P.6-3 |
| 2   | "6.1.2 Action Bar"             |  P.6-4 |
| 3   | "6.1.3 Search Criteria"        |  P.6-5 |
| 4   | "6.1.4 Study List"             |  P.6-5 |
| 5   | "6.1.5 Series Tab"             |  P.6-7 |
| 6   | "6.1.6 Options Tab"            |  P.6-7 |
| 7   | "6.1.7 Messages & Information" |  P.6-8 |

### 6.1.1 Study Location

The [Location] box and buttons [System] and [Drives] set the target of the [Search] button as follows:

| No. | Display   | Purpose  |
|-----|---|--|
| 1   |   | [System] sets the [Search] target to the local system hard drives plus all configured servers. Enter search criteria and click [Search] to find studies. |
| 2   |  | [Drives] sets the [Search] target to the CD/DVD drives. It then searches all CD/DVD drives, displaying the results in the Study List.                    |
| 3   |  | [Location] becomes a drop-down list enabling individual selection of servers and drives as the [Search] target.  |










## 6.1.2 Action Bar



Studies Management functionality is accessed via the action bar.



Not all Action Bar buttons appear at all times and in all modes.


The purpose of each Action Bar button is summarized as follows:

| No. | Button  | Name        | Purpose   |
|-----|---|-------------|---|
| 1   |    | Search      | Searches based on the search target set by the [System] and [Drives] buttons, plus all search criteria. All matching studies are displayed in the Study List. The Date search criteria is always used.  |
| 2   |    | View        | Displays the selected study or studies in the Image Viewer window. Alternatively, you can double-click anywhere in a Study List row to view a study.  |
| 3   |   | Dose Report | (Available if there is dose report.)<br>Displays the dose report.   |
| 4   |  | Series      | (For non-local studies) When a single study is selected in the Study List, Series displays all study series types contained in this study. When two or more studies are selected in the Study List, this button causes the Study Series information to be retrieved for each study but not displayed until only a single study is selected. Series are automatically displayed for local studies.   |
| 5   |  | Receive     | Available only for studies on a network server or CD/DVD, [Receive] copies all selected studies to the local system.  |
| 6   |  | Send        | Sends selected studies to another server or system. If configured, notification messages are displayed and emails are sent for each send.   |
| 7   |  | Information | Available only when a single study is selected in the Studies Management window for the local system or CD, [Info] displays detailed information about the patient who is the subject of the selected study. You can also use this button to anonymize or modify information in the local study.  |
| 8   |  | Delete      | (Present for users only if enabled by installation personnel) Available only for unprotected studies on the local system, [Delete] permanently deletes the selected studies from the local system after asking for confirmation. Protected studies cannot be deleted, although they can be first unprotected (see [Unprotect] below) and then deleted. To preserve the integrity of your archive system, do not delete studies that have not yet been archived to server or CD/DVD. |
| 9   |  | Write       | Available only for studies on the local system or a network server, [Write] prepares the selected studies for writing and then writes them to CD/DVD.   |

| No. | Button  | Name      | Purpose   |
|-----|---|-----------|---|
| 10  |  | Protect   | Available only for studies on the local system, [Protect] enables protection for the selected studies so that they cannot be deleted.   |
| 11  |  | Unprotect | Available only for studies on the local system, [Unprotect] clears protection for the selected studies. It is generally not recommended to unprotect studies that have not yet been archived to server or CD/DVD. |













Detailed usage information for these buttons is provided as needed throughout this manual.

### 6.1.3 Search Criteria

One Search Criteria box is provided above each column in the Study List. You can filter what is searched for by entering words or parts of words in one or more Search Criteria boxes and then clicking [Search]. All studies matching the search criteria are listed. See  ["6.2.2 Finding Studies" P.6-10.](#)

### 6.1.4 Study List

The Study List provides one row of information per study in columns organized as follow:

| No.   | Column                        | Contents  |      |                               |   |                           |   |                             |   |                     |   |             |
|---|-------------------------------|---|------|-------------------------------|---|---------------------------|---|-----------------------------|---|---------------------|---|-------------|
| 1   | Icon                          | <div>Denotes each study's location and protection state as follows:</div> <table><thead><tr><th>Icon</th><th>Location and Protection State</th></tr></thead><tbody><tr><td></td><td>Protected on local system</td></tr><tr><td></td><td>Unprotected on local system</td></tr><tr><td></td><td>On a network server</td></tr><tr><td></td><td>On a CD/DVD</td></tr></tbody></table> | Icon | Location and Protection State |  | Protected on local system |  | Unprotected on local system |  | On a network server |  | On a CD/DVD |
| Icon  | Location and Protection State |   |      |                               |   |                           |   |                             |   |                     |   |             |
|  | Protected on local system     |   |      |                               |   |                           |   |                             |   |                     |   |             |
|  | Unprotected on local system   |   |      |                               |   |                           |   |                             |   |                     |   |             |
|  | On a network server           |   |      |                               |   |                           |   |                             |   |                     |   |             |
|  | On a CD/DVD                   |   |      |                               |   |                           |   |                             |   |                     |   |             |
| 2   | Accession Number              | A unique identification number for the study.   |      |                               |   |                           |   |                             |   |                     |   |             |
| 3   | Patient Name                  | The full patient name, with last name first.  |      |                               |   |                           |   |                             |   |                     |   |             |
| 4   | Patient ID                    | The patient ID code.  |      |                               |   |                           |   |                             |   |                     |   |             |
| 5   | Physician Name                | The full name of the study's performing physician, with last name first.  |      |                               |   |                           |   |                             |   |                     |   |             |
| 6   | Date                          | The date on which the study was performed.  |      |                               |   |                           |   |                             |   |                     |   |             |
| 7   | Origin                        | An institution-defined field which more precisely indicates the origin of the study. For example, this could indicate the diagnostic lab in which the study was performed. There can be multiple Origins per Site.  |      |                               |   |                           |   |                             |   |                     |   |             |
| 8   | Modality                      | The DICOM abbreviation for the imaging modality, for example, XA for angiography.   |      |                               |   |                           |   |                             |   |                     |   |             |

### ■ Sorting

Click any column heading to sort by information in the column. Each click of a column heading alternates between ascending and descending sort order.

### ■ Scrolling in the Study List

You can scroll the Study List contents both horizontally and vertically by clicking and dragging the respective scroll bars.

The vertical scroll bar only appears when there are too many studies to fit in the Study List display area. To scroll the Study List vertically, drag (point to slider, click down without releasing, and drag mouse in desired direction) the vertical scroll bar slider downward to see studies further down the list, or upward to see studies up the list. Alternatively, spin the mouse wheel (if available) toward yourself to scroll down the list or away from yourself to scroll up the list.

The horizontal scroll bar only appears when all columns of information cannot fit across your display. To scroll the Study List horizontally, drag the horizontal scroll bar slider to the right to see more columns to the right or drag it to the left to see more columns to the left. If you have a high-resolution display, the horizontal scroll bar may not appear at all.

### ■ Widen or Narrow Columns as Desired

If the text in any columns ends with "..." that means there is hidden text to the right. To reveal the text, use the mouse to widen the column by pointing to its right heading edge, waiting for the mouse pointer to change to a vertical bar with arrows, and then clicking the mouse and dragging the column edge to the right. Release the mouse button once the column is the desired width. To narrow a column, do the opposite by dragging the heading's right column edge to the left.

| Patient ID        | Physician |
|-------------------|-----------|
| PIDDAR9500REV_011 |           |
| PIDDAR9500REV_... |           |


### ■ Making Selections in the Study List

You can select studies in the Study List in a variety of ways as follows.

- To select a single study, click the mouse anywhere in the desired study row. The row becomes highlighted.
- To select several non-adjacent studies, use this [Ctrl]-click technique. Point to the first study row and click the left mouse button. The row becomes highlighted. Press and hold down the [Ctrl] key and then continue clicking all other desired rows without releasing the [Ctrl] key. When finished selecting rows, release the [Ctrl] key. All selected rows are highlighted.
- To select two or more adjacent study rows, use this [Shift]-click technique. Point to the highest row that you want and click the left mouse button. Point to the lowest row that you want and hold down the [Shift] key and click the left mouse button. Release the [Shift] key. All selected rows are highlighted.
- To select all studies in the Study List, use [Ctrl]-A (press [Ctrl] key and [A] key simultaneously).




### 6.1.5 Series Tab

Sharing the same screen area as [Options] tab, the Study Series appears when you click the [Series] tab to the left of the [Options] tab. The [Series] tab is shown by default. Most studies contain one or more series of images, image loops or special items such as reference images. When a single local study is selected in the Study List, all series for the study are shown in the Study Series with one icon per type. Series type icons include a brief label beneath them, indicating the type preceded by the number of items within that series.



If MPPS support is available, specific series detail can be displayed. Select one series and click [Series Detail] at the right edge of [Series] tab.

Possible series types and their corresponding icons are as follows:

| No. | Icon  | Series Type                      |
|-----|---|----------------------------------|
| 1   |   | Cine (raw)<br>Cine DSA (raw)     |
| 2   |  | Photo Capture<br>Reference Image |
| 3   |  | Dose Report<br>Dose Report Image |



The Studies Series type icons are only displayed when a single study is selected in the Study List.


### 6.1.6 Options Tab

Sharing the same screen area as Study Series, the [Options] tab appears when you click its title to the right of the [Series] tab in the lower-left window area.


 ["17.7.2 Options Tab"](#)

### 6.1.7 Messages & Information


The Messages & Information area at the bottom of the window provides messages, error messages, DICOM UID (Unique Identifier), study count, a message log indicator, a network activity indicator, and software version information. For detailed network activity information, click the [Network Activity] indicator, an information box appears. Click the [Network Activity] indicator again to close the information box.



For details of version license information, click [Software Information]. Information dialog box of software is displayed.



To display the license information, click [License Info] button. The license information dialog box is displayed.













## 6.2 Managing Studies

### 6.2.1 Introduction

This section describes how to perform common study management tasks as follows.

| No. | Section                                  | Page   |
|-----|--|--|
| 1   | "6.2.2 Finding Studies"                  |  P.6-10 |
| 2   | "6.2.3 Downloading Studies from Servers" |  P.6-12 |
| 3   | "6.2.4 Displaying RDSR"                  |  P.6-13 |
| 4   | "6.2.5 Importing Studies from CD/DVD"    |  P.6-14 |
| 5   | "6.2.6 Sending Studies"                  |  P.6-14 |
| 6   | "6.2.7 Writing to CD/DVD"                |  P.6-15 |
| 7   | "6.2.8 Setting Study Protection"         |  P.6-19 |
| 8   | "6.2.9 Deleting Local Studies"           |  P.6-20 |



Before continuing, become fully familiar with  "6.1 The Studies Management Window" P.6-2.

## 6.2.2 Finding Studies




If you open the Studies Management window while a new study is in acquisition, a find is automatically performed for all studies with the same Patient ID as the active study.

To find studies to work with, follow this procedure from the Studies Management window.

- 1 Set the Date criteria to include the dates of the desired studies. Choose the desired value from the Date list.


| No. | Date         | Selection   |
|-----|--------------|---|
| 1   | Today        | Studies dated today.  |
| 2   | Yesterday    | Studies dated yesterday.  |
| 3   | This Week    | Studies dated in this calendar week.<br>For example, if you choose this on a Wednesday, the Studies list will be restricted to the four days Sunday through Wednesday.  |
| 4   | Last 7 Days  | Studies dated within the last 7 days including today.   |
| 5   | This Month   | Studies dated on any day so far in this calendar month, including today.<br>For example, if you choose this on July 11, the Studies List will be restricted to studies dated July 1 through 11.   |
| 6   | Last 30 Days | Studies dated within the last 30 days, including today.   |
| 7   | This Year    | Studies dated on any day so far in this calendar year, including today.<br>For example, if you choose this on January 15th, the Studies List will be restricted to 15 days, whereas if you choose this on December 15, the Study List will be restricted to approximately 349 days. |
| 8   | Last Year    | Studies dated on any day in the previous calendar year.   |
| 9   | All Time     | All dates.<br>When a server is included in the search, it is enforced that the All Time Date criteria be combined with at least one other criteria item to avoid excessive results.   |
| 10  | Custom Date  | Studies dated on the specific date chosen from a calendar.  |

- 2 If desired, click the arrow to the left of Accession Number to show additional search criteria fields. This is available only when MPPS is enabled.  ["DICOM Tab" P.17-24.](#)

|                     |                     |                    |                         |                           |                     |          |
|---------------------|---------------------|--------------------|-------------------------|---------------------------|---------------------|----------|
| Accession Number    | Patient Name        | Patient ID         | Physician Name          | Date                      | Origin              | Modality |
|                     |                     |                    |                         | 02/20/2018                |                     |          |
| Study Time          | Study ID            | Study Instance UID | Modalities in Study     | Referring Physician Names |                     |          |
| 00:00:00 - 23:59:59 |                     |                    |                         |                           |                     |          |
| Series Number       | Series Instance UID | Request Proc. ID   | Scheduled Proc. Step ID | Performed Proc. Step Date |                     |          |
|                     |                     |                    |                         | 02/20/2018                | 00:00:00 - 23:59:59 |          |
| Instance Number     | SOP Instance UID    |                    |                         |                           |                     |          |
|                     |                     |                    |                         |                           |                     |          |

- 3 Enter your search by entering whole or partial words that you know appear in the desired studies in one or more search criteria boxes.

By default, search criteria matching is only performed at the beginning of the searched text in the corresponding fields of the Search List (for example: criteria\*). This can be changed so that search criteria matching is performed throughout the searched text (for example:

\*criteria\*) by checking the [Search by keyword] option as described in  ["17.7.2 Options Tab" P.17-14.](#)

- 4 Click [Search] to perform the search.


The Studies List is now restricted to all studies that match the search criteria.


**Studies Management**

Location: System | System | Drives

Search | View | Series | Send | Info | Delete | Write | Protect | Unprotect | Close

| Accession Number | Patient Name  | Patient ID         | Physician Name | Date                     | Origin        | Modality |
|------------------|---------------|--------------------|----------------|--------------------------|---------------|----------|
|                  |               |                    |                | 02/18/2018 to 02/22/2018 |               |          |
| Accession Number | Patient Name  | Patient ID         | Physician Name | Date                     | Origin        | Modality |
|                  | Patient Cindy | PIDDAR9500REV_011  |                | 02/22/2018 15:41:52      | DAR9500REV-01 | XA       |
|                  | Patient Bryan | PIDDAR9500REV_0... |                | 02/22/2018 15:41:08      | DAR9500REV-01 | XA       |
|                  | Patient Anne  | PIDDAR9500REV_0... |                | 02/22/2018 15:34:10      | DAR9500REV-01 | XA       |

If there are too many studies to fit in the Study List, use the vertical scroll bar to scroll downward as previously described in  ["Scrolling in the Study List" P.6-6.](#)

- 5 To see columns that have scrolled off the right side of the Study List, use the horizontal scroll bar.
- 6 If a column contains text that ends in "..." widen the column as previously described in  ["Widen or Narrow Columns as Desired" P.6-6.](#) If this causes a column to disappear off the right edge, use the horizontal scroll bar to scroll to the right.

### 6.2.3 Downloading Studies from Servers


Once you have found the desired studies in the System Study List, you can download them to the local system.




#### NOTE

If you do not download studies before writing them to CD/DVD, they will first be automatically downloaded to the local system before any writing occurs. You may find it preferable to manually download the studies as described here to avoid delaying the writing process.

To download studies from the network, follow this procedure.

- 1 Click [System] button.
- 2 Enter the desired search criteria and click [Search].  
All matching studies contained on the local system and all configured servers are displayed. Studies on a server are denoted with  in the Icon column of the Study List.
- 3 With the desired studies selected in the Studies List, click the [Receive] button. A confirmation box appears.
- 4 Click [OK] to confirm the downloading of all highlighted files.



- 5 A message is displayed once a study has been downloaded.

### 6.2.4 Displaying RDSR

If there is RDSR in the study, click [RDSR] on the action bar to display the [Dose Report] window to check the description.

| Dose Report |        |          |                            |             |          |              |         |           |    |
|-------------|--------|----------|----------------------------|-------------|----------|--------------|---------|-----------|----|
| Acq.        | Frames | Time     | DUP                        | DAP (mGym2) | AK (mGy) | Plane        | LAO/RAO | CAUD/CRAN | kV |
| 1           | 9      | 13:50:08 | Fluoro_Cardio_15pps        | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
| 2           | 11     | 13:50:15 | Rad_Cardio_CAG[15f-10s]    | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
| 3           | 15     | 13:50:23 | Rad_Cardio_DSA[12f-20s]    | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
| 4           | 12     | 13:50:33 | Rad_Peri_DSA[3f-2f-1f]     | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
| 5           | 25     | 13:50:42 | Rad_Peri_DSA-HS-FlexAPS    | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
| 6           | 27     | 13:50:51 | Rad_Peri_RSM[12f-45s]-S... | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
| 7           | 17     | 13:51:08 | Rad_Peri_RSM[12f-45s]-S... | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
| 8           | 26     | 13:51:18 | Rad_Cardio_CAG[15f-10s]    | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
| 9           | 21     | 13:51:26 | Rad_Cardio_CAG[15f-10s]... | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
| 10          | 23     | 13:51:33 | Rad_Cardio_CAG[15f-10s]... | 0.000       | 0.000    | Single Plane | LAO 0   | CRAN 0    | 40 |
|             |        |          |                            |             |          |              |         |           |    |
|             |        |          |                            |             |          |              |         |           |    |
|             |        |          |                            |             |          |              |         |           |    |

Cumulative DAP (mGym2): n/a  
Cumulative AK (mGy): n/a  
Fluo Time (min.): 0.0

Plane: ☐ Frontal ☐ Lateral ☐ Both  
Type: ☐ Fluoroscopy ☐ Radiography ☐ Both  
☐ Display Pulsed Fluoroscopy Information as Continuous Fluoroscopy

Restore  
Close

On the [Dose Report] window, enable to modify the display by following operations.

- Drag and move the title to change the order of titles.
- Right-click on the title to display the list to select the display title.
- When check [Display Pulsed Fluoroscopy Information as Continuous Fluoroscopy], enable to display "mA", "ms" and "Exposure time" as continuous fluoroscopy.
- Select desired plane in [Plane] to display dose information of the plane.
- Select Fluoroscopy or Radiography in [Type] to display dose information of the type.
- Click [Restore] button to return to the initial setting.
- Click [Close] button to save the changes.



#### NOTE

Method of calculation for fluoroscopy time is different on RDSR and system. Therefore, fluoroscopy time is not correspond completely with each other.

For RDSR, calculate with pulse length and pulse number to consistent between each parameter. By contrast, time the fluoroscopy time at regular time interval between turn on and off the foot switch of high-voltage generator.

### 6.2.5 Importing Studies from CD/DVD

You can import studies to the local system from any compatible DICOM Standard CD/DVD. This often useful when you wish to work with studies that were performed outside your institution.


To import studies from CD/DVD, follow this procedure.


- 1 Insert the CD/DVD from which you wish to import studies into the CD/DVD drive.
- 2 Click the [Drives] button in the upper-left corner of the Studies Management window.
- 3 The CD/DVD drive is searched for DICOM studies. All studies found appear in the Study List of the Studies Management window.  
  
A warning message is displayed if the CD/DVD drive contain no DICOM studies. This message can be ignored if the drive is empty or contains a disc.
- 4 To import one or more studies from CD/DVD into the local system.  
Select them in the Drive Study List and then click [Receive]. A confirmation box appears.
- 5 Click [OK] to confirm the importation of all highlighted files. Once completed, a message to that effect is displayed.
- 6 You can now see the newly-imported studies in the Studies Management window. Click [System] to return to the Studies Management window. Adjust or remove any search criteria to that the studies can be seen.


### 6.2.6 Sending Studies

Studies, loops, images and dose report can be transmitted to DICOM station. To sent studies, follow this procedure.

- 1 Select the desired studies in the Study List.
- 2 Click [Send]. A menu appears with a choice of destinations.



- 3 Select the sent type and click the desired destination.  
Check [Process Before Transfer] to start image processing on playback frame and transfer.  
See  "DICOM Hosts Tab" P.17-20 for processing setting.
- 4 When the confirmation prompt appears, click [OK] to confirm the operation.

If there is no [Send] button on the top action button bar, a destination host must first be configured as described in the DICOM Hosts AE Table item of the  "DICOM Hosts Tab" P.17-20.

If configured, an administrator-defined message will be displayed on the local system upon completion of each send operation. One message appears per send operation whether the operation includes one study or may.

Once all studies that you wish to write are visible in the System Study List, you can select them and then write them to CD/DVD.

### 6.2.7 Writing to CD/DVD

#### Recommended Media (As of January, 2015)

\*Recommended media may change in the future without notice.

Disc, CDR80WKY10SV (Parts No.: 088-58818-18)

Disc, DR-47WWY10SNT (Parts No.: 088-58818-21)



#### CAUTION

Use the media which is provided from Shimadzu for the system. There is no warranty for functionality and operation when using other media.



#### NOTE

Unable to use DVD-RW.



#### NOTE

Do not write to the media during a study. Write to the media after completing the study.



Recommended media as of January 2012. The recommended Media may change in the future without notice.



Writing may stop when start writing on DVD after writing on CD.


### Select the Studies

- 1 Select desired studies as described in ["6.1.4 Study List" P.6-5](#) or ["" P.6-7](#).
- 2 Click the [Write] button.  
You are prompted to insert blank disk into the writer. Do this and close the disc tray. Any needed file conversion is performed.



If the [Auto-accept Default] option is configured (see ["Devices Tab" P.17-56](#)), writing begins and you are only prompted for additional blank discs. Continue with section ["Insert Additional Blank CD/DVDs" P.6-19](#).

- 3 When the Media Writing dialog box appears, optionally change the proposed volume name. You can select the writing item from [Images, [RDSR] or [Both]. Enter a volume name and select the writing item and then click [Start] button.  
The proposed name is based on a unique sequence number so if you change the volume name, it is advisable to keep the volume name unique, possibly by including numerical digits.






## Optionally Select Fewer Items

If you wish to write fewer items than selected earlier on the Studies Management window, follow this procedure.

### 1 Click [Select File].

The detailed selection window appears, initially showing one line of information per study with three columns of information: Accession Number, Patient Name and Patient ID. This information is displayed in a tree format with the first level being the study, the second level being the study series type, and the third and final level being the study series items themselves.




In this example, 3 studies are shown. Any tree branch can be expanded (click the "+" symbol) to reveal all series types for the study, and each series type branch can be expanded to reveal the items in the series.

### 2 Depending on what level of the tree you click, you can select or deselect studies (first level), study series types (second level), or the actual series items themselves (third level).


For the first two levels, a "+" or "-" symbol appears at the left edge of each line in the tree. Click the "+" symbol to expand that branch. Click the "-" symbol to contract an expanded branch.

When the tree-based selection window first appears, every study is selected as indicated by the gray highlight. Click "+" symbols as needed to reveal study series types and the series items. Click "-" symbols to contract studies or series types you do not need to see in detail. Once you start clicking the "+" and "-" symbols, the global selection disappears and the individual studies, study series types, and series items for the line in the tree that you clicked.



- 3 You can now use the [Shift]-click techniques (selecting adjacent item) or [Ctrl]-click techniques (selecting non-adjacent item) to select specific studies, study series types, and series items.

In this example, one of the studies is only partially selected.



- 4 Observe the Space used percentage valued at the bottom of the selection window. Percent values greater than 100 indicated that two or more discs must be written. For space efficiency, it is recommended that you make selections so that the percent value is close to but not over 100 % or 200 %, and so on.

### Perform the Write

Follow this procedure to complete the writing process.

- 1 To begin writing, click the [Start] button. All study series items are checked for integrity.  
  
If any errors are detected and the Auto-accept Defaults option is not configured, the study series item tree window of the Media Writing dialog box appears with the offending items highlighted. If desired, take note of the problem items. Click [Continue] to go ahead and write everything else except the offending study items.  
  
If errors are detected and the Auto-accept Defaults option is configured, automatic continuation occurs without prompting, skipping items with errors.  
  
CD/DVD writing begins. During the write process, the CD/DVD writer tray is locked shut. Do not attempt to open the tray until the CD/DVD writing finishes.
- 2 If a defective blank disc or a non-blank disc is inserted, it will be rejected and you will be prompted to insert a new blank disc. Discard the rejected disc and insert a new one.  
  
Close the disc tray to continue.



By default, DSA (Digital Subtraction Angiography) images are saved to media in their subtracted state. This enables easy DSA image viewing in most viewers.

### ■ Insert Additional Blank CD/DVDs

When CD/DVD writing finishes, you may be prompted to insert another blank disc. Keep inserting blank discs when prompted until finished.

When writing is completed, insert the media to tray again. Display from the Studies Management window and check that CD/DVD has written correctly.

### ■ Cancel a Media Write Operation

Usually, if you cancel a media (CD/DVD) write, the media must be discarded. To cancel a media write, click [Cancel Write] on the Image Viewer window on the side menu. Or, in the Studies Management window, click [Write] and then respond [Yes] to the [Cancel Media Write] prompt.

### ■ Label CD/DVD

According to your organization's standard, label the CD/DVD, typically at least including the volume name and date. CD/DVD-safe permanent markers can be used to write directly on the top of the CD/DVD.



If you plan on affixing paper labels to the CD/DVD, be aware that the labels could partially come off while inside a writer, possibly causing damage, or the CD/DVD may vibrate when spinning due to an unbalanced label. If you must use such labels, assure that they are perfectly centered and are firmly attached without any loose edges.

See "16.1 The Mini Viewer" P.16-2 for mini viewer.

## 6.2.8 Setting Study Protection

Studies can be protected or unprotected as desired by respectively clicking the [Protect] or [Unprotect] action buttons. The first column of the Study List indicates protection status, either protected (closed-lock icon ) or unprotected (open-lock icon ).

### 6.2.9 Deleting Local Studies

---


(Feature present for uses only if enabled by installation personnel.) Studies on the local system that are in an unprotected state can be manually deleted with the [Delete] button. To manually delete a study, select one or more unprotected studies and click [Delete]. Confirm or cancel the deletion.



Studies on a network server cannot be deleted. Protected studies cannot be deleted until they have been unprotected. Eventually, unprotected studies that are not manually deleted will be automatically deleted if the system needs to recover the cache space that they are occupying. Automatic deletion of unprotected studies occurs in least-recently accessed order.

# Chapter 7

## Review Reference

Become fully familiar with  ["4.2 The Main User Interface Window"](#) [P.4-2](#) before beginning this chapter.

For as long as a study remains on the local system (including studies viewed from the network or CD/DVD), the following settings are preserved on a per loop/image basis: Brightness/Contrast, Auto Window Level, LUT, Filters, Zoom, Pan, and DSA parameters such as shutter positions and mask number.

This chapter provides detailed reference information on the Review process as follows.

### Description

|                      |  |                      |
|----------------------|--|----------------------|
| <a href="#">7.1</a>  | <a href="#">Controlling Image Playback . . . . .</a>     | <a href="#">7-2</a>  |
| <a href="#">7.2</a>  | <a href="#">Adjusting Image Appearance . . . . .</a>     | <a href="#">7-4</a>  |
| <a href="#">7.3</a>  | <a href="#">Working with Multiple Studies . . . . .</a>  | <a href="#">7-9</a>  |
| <a href="#">7.4</a>  | <a href="#">Displaying Study Information . . . . .</a>   | <a href="#">7-10</a> |
| <a href="#">7.5</a>  | <a href="#">Modifying Study Information . . . . .</a>    | <a href="#">7-12</a> |
| <a href="#">7.6</a>  | <a href="#">Anonymizing Study Information . . . . .</a>  | <a href="#">7-15</a> |
| <a href="#">7.7</a>  | <a href="#">Working with Image Annotations . . . . .</a> | <a href="#">7-17</a> |
| <a href="#">7.8</a>  | <a href="#">Saving Images . . . . .</a>                  | <a href="#">7-20</a> |
| <a href="#">7.9</a>  | <a href="#">Print Images . . . . .</a>                   | <a href="#">7-21</a> |
| <a href="#">7.10</a> | <a href="#">Modifying Image Comments . . . . .</a>       | <a href="#">7-33</a> |

# 7.1 Controlling Image Playback

As soon as a study is opened for review, the first loop begins playback. You control playback via the mouse, keyboard or IVR NEO/IVR Shuttle/SMART Touch.

## 7.1.1 Mouse and Keyboard

The Cine Control bar on the bottom border of the Image Viewer window provides seven GUI buttons for controlling playback with mouse clicks.



The keyboard arrow keys and mouse wheel can also be used.

| No. | Action   | Keyboard | GUI |
|-----|--|----------|-----|
| 1   | Play/show next loop/image.   |          |     |
| 2   | Play/show previous loop/image.   |          |     |
| 3   | If playing, pause.   | or       |     |
| 4   | If paused, show next/previous frame.   | /        | /   |
| 5   | If paused, start playing. <i>(Keyboard only: Hold key for more than 2 seconds to begin forward / backward play.)</i> | /        | /   |



Use only the arrow keys on the keyboard edit group to the right of the main keys.



Although intended for use during acquisition, the IVR NEO/IVR Shuttle/SMART Touch is capable of many playback functions.


["4.8 IVR NEO \(Either-or SMART Touch\)" P.4-65](#), ["4.10 IVR Shuttle \(Option\) \(Either-or SMART Touch\)" P.4-73](#), ["4.11 SMART Touch \(Either-or IVR NEO/IVR Shuttle\)" P.4-77](#)

### Play Acquisition Monitor with Keyboard

Hold down [Ctrl] key on keyboard and use arrows keys above to control the play of acquisition monitor. Release [Ctrl] key to switch the play to reference monitor.


7.1.2 Loop Playback Mode




Use the button on the side menu to change loop playback mode.



Change the mode "Play Single Loop" and "Play Single Loop with Synchronize".

Three Loop Playback Mode buttons are available in the bottom bar.



| Button  | Name and Function  |
|---|--|
|  | <b>Play All Loops:</b> Plays back each loop, one after the other. Still images are skipped. Repeats once all loops have been played.   |
|  | <b>Play Single Loop:</b> Plays back the selected loop continuously.  |
|  | <b>Play Single Loop with Synchronize:</b> <i>(Reference monitor only)</i><br>Like Play Single Loop, plays back selected single loop until a new loop is received. Then, the new loop is selected and played back continuously. |

7.1.3 Adjusting Playback Speed

The current playback speed, in frames (images) per seconds (fps) is displayed at the bottom of the Play Speed button on the side menu. This button enables you to precisely adjust loop playback speed in either direction (from 10 % to 200 % of recorded speed) as follows.



- Right-click the Play Speed button once to increase the playback speed slightly. Right-click it multiple times to increase playback speed gradually up to maximum speed. Right-click and hold the mouse button to rapidly increase playback speed.
- Conversely, left-click the Play Speed button to reduce playback speed.
- To reset playback speed to recorded speed, click the Play Speed button simultaneously with both mouse buttons.

## 7.2 Adjusting Image Appearance

You adjust image appearance using GUI buttons and the keyboard.




The monitors are already pre-adjusted to optimize contrast and brightness. Do not adjust these settings on the monitors. Instead, make your adjustments as described here. Image-appearance adjustments do not alter image files, although the effects of the adjustments are preserved per image/loop for as long as the study is on the local system.

### 7.2.1 Basic Brightness/Contrast Control

Image brightness and contrast are adjusted on screen using the respective side menu buttons or with the mouse.

#### ■ Using Side Menu Buttons



Adjust image brightness (left button) as follows:

- Right-click the Adjust Brightness button to increase image brightness (maximum 100).
- Conversely, left-click the Adjust Brightness button to decrease image brightness (minimum -100).
- To reset brightness to its default, click the Adjust Brightness button simultaneously with both mouse buttons.

Adjust image contrast (right button) as follows:

- Right-click the Adjust Contrast button to increase image contrast (maximum 2000).
- Conversely, left-click the Adjust Contrast button to decrease image contrast (minimum -100).
- To reset contrast to its default, click the Adjust Contrast button simultaneously with both mouse buttons.

#### ■ Using the Mouse

Change image brightness and contrast as follows:

- Press and hold both left and right buttons.
- Move the mouse up/down (brightness) or left/right (contrast) within the visualization zone.


The change in image brightness and contrast is proportional with the mouse speed. Keeping the [Alt] key pressed makes the change independent of the mouse speed.



## 7.2.2 Advanced Brightness/Contrast Control

### ■ Auto Window Level Feature

When enabled, click [Filters] button on the side menu and select [Auto W] button from the fly-out menu. The Auto Window Level Feature enables brightness/contrast optimization on a frame-by-frame basis.



This automatically compensates for image loops whose images have varying brightness and contrast levels from frame to frame. This can occur, for example, when some frames from a loop of a beating heart, feature many bright areas whereas other frames are substantially darker. When Auto Window Level is enabled, the Brightness and Contrast controls are used to fine tune the brightness/contrast (-100 to 100 Brightness, -100 to 2000 Contrast) set by Auto Window Level. For example, if the results produced by Auto Window Level appear too dim, increase the brightness by right-clicking the brightness control.



There are two sets of Brightness and Contrast values; one for when Auto Window Level is active, and one for when it is not.


## LUT Feature

To further optimize image brightness and contrast, you can apply special brightness/contrast profiles to your images via the LUT (Look-up Table) feature.

To apply a LUT profile to our images, follow this procedure.

- 1 Click [Filters] on the side menu.

The LUT fly-out menu appears.



- 2 Click one of the Curve buttons [LUT-A], [LUT-B], or [LUT-C] to select one of three standard profiles.

- LUT-A: Darkens the image by making all grays except for those close to black or white, lighter.
- LUT-B: Darkens the image by making all grays except for those close to black or white, darker.
- LUT-C: Increases the contrast by darkening the dark grays and lightening the light grays while having little effect on the middle grays.

- 3 Alternatively, click [LUT-O] and choose the desired custom profile from the list.

A graph above the list indicates approximate profile characteristics.

- 4 To revert to the default LUT profile, click [Linear] on the LUT fly-out menu.

Additional LUT profiles can be created by the administrator as described in [👉 "17.5 Working with LUT Profiles" P.17-4.](#)




### NOTE

The chosen LUT profile applies only to the selected loop/image.

### 7.2.3 Applying Filters

Image appearance can be enhanced by applying filters. Click [Filters] on the side menu. Several buttons fly out to the right. Click one of the first three buttons to set a filter. A filter is only active when its button appears pushed in. Click [Reset] to deactivate all filters. The filters are used as follows.



- The first two filters, Sharpness (Unsharp Mask) and Convolution (High Pass Convolution) enable you to increase the perceived sharpness of the image through edge enhancement. Increase the effect by right-clicking the button or decrease the button in the range -100 to 100. Only one of these filters can be active at a time. Click the button simultaneously with both mouse buttons to set the filter to its default (fixed at 0).
- Click [Negative] to invert the image so that black is changed to white, white is changed to black, dark gray changed to light gray, and so on. This effect is in addition to either sharpness filter. To clear this filter, click the button again so that it no longer appears pushed in.
- Images can be flipped horizontally or vertically. Click the [H Flip] button (the figure on the left) or [V Flip] button (the figure on the right) on the side menu to flip horizontally or vertically.



### 7.2.4 Zooming and Panning

---

You can zoom and pan images with the IVR NEO/SMART Touch, the GUI and mouse, or the keyboard.

#### Using the Zoom Button and Mouse

The current zoom level is displayed at the bottom of the [Adjust Zoom] button on the side menu.



This button enables you to precisely control zoom level as follows.


- Right-click the Adjust Zoom once to increase the zoom level slightly. Right-click it multiple times to increase zoom level gradually up a maximum of 2.5 (250 %) for 1024x1024 images and 5 (500 %) for 512x512 images. Right-click and hold the mouse button to rapidly increase zoom level.
- Conversely, left-click the [Adjust Zoom] button to reduce zoom level. Minimum zoom is 1 (100 %) for 1024x1024 images and 2 (200 %) for 512x512 images.
- When zoomed in, you can pan the image by pointing to the area of interest, clicking down on the right mouse button, and without releasing the button, moving the mouse pointer toward screen center. Release the mouse button once the area of interest is centered on screen.
- To reset zoom level to the default (1 for 1024x1024 images, and 2 for 512x512 images), click the [Adjust Zoom] button simultaneously with both mouse buttons.

## 7.3 Working with Multiple Studies


It can be convenient to work with more than one study at a time.

Follow this procedure.


- 1 In the Image Viewer window, click [Close] to close any studies and then click [Open] to open the [Studies Management] window.
- 2 Select multiple studies in the [Studies Management] window.  
Hold down [Shift] key to select on the range and [Ctrl] key to select individually.




- 3 Click [View].  
The newest study's first loop is played in the Image Viewer. Also a [Selected Studies] list is displayed at the bottom of the window with a tick mark next to the selected study.



- 4 Choose loops and images for display as desired.
- 5 Click a different study in the [Selected Studies List].  
A tick mark appears next to the study and the study's first loop is played in the Image Viewer.



- 6 Display loops/images and select other studies as desired.
- 7 Click the [Hide] button to hide the Selected Studies list.  
Click the [Selected Studies List] icon  on the bottom window border to re-display it.

## 7.4 Displaying Study Information

As seen in the chapter "5 Review Basics", a [Study Information] dialog box can be displayed by clicking the [Info] button in the Studies Management window. The same information window can be displayed by clicking [Study] on the side menu, or by clicking the [Display Study Information] button on the right side of the bottom window border and then clicking the [Display Study Information] button (top) on the fly-put menu that appears.



The [Study Information] dialog box includes patient demographics and other information:

**Study Information**

Patient

Patient ID: P1DDAR9500REV-022 Date of Birth: 07/07/1940 Age: 79

Patient Name: Prefix: John Middle Name: Last Name: Shimadzu Suffix:

Height: inches Weight: lbs Sex: Male Select Orientation

Study: 07/12/2019 14:18:28

Study ID: S00001 Accession Number: A00001

Modality: LA SR Referring Physician: Referring First

Study Description:

Performing Physician's List: [x] Performing First Operators' List: [x] Operator Second

Primary Physician: Performing First

Anonymize Modify Close

Expert users can also display DICOM header information by clicking the Display Study Information button, and then clicking the Display all DICOM Information button (bottom) on the fly-out menu that appears.





This information is mainly of interest to those familiar with interpreting DICOM image file headers.

The DICOM Information window provides detailed DICOM header information:



## 7.5 Modifying Study Information

Text-based study information, including patient demographics, can be modified locally for the purpose of making minor corrections before printing or exporting images.



Modifications are only made on the local system.

To modify study information, follow this procedure.

- 1 Click [Open] button on the side menu to pen the Studies Management window.
- 2 Click the study you wish to modify.
- 3 Click [Info]. The Study Information dialog box appears.
- 4 Click [Modify] and enter the administrator (super) password, and click [OK]. The [Study Information Modification] dialog box appears.

| Study Information Modification  |                             |                                  |          |      |
|---------------------------------|-----------------------------|----------------------------------|----------|------|
| Patient Information             |                             | Patient Name                     |          |      |
| Patient ID                      | PIDDAR9500REV_022           | Prefix                           |          |      |
| Height                          | <input type="text"/> inches | First Name                       | John     |      |
| Weight                          | <input type="text"/> lbs    | Middle Name                      |          |      |
| Sex                             | Male                        | Last Name                        | Shimadzu |      |
|                                 |                             | Suffix                           |          |      |
|                                 |                             | Date of Birth                    |          |      |
|                                 |                             | Year                             | 1940     |      |
|                                 |                             | Month                            | 07       |      |
|                                 |                             | Day                              | 07       |      |
|                                 |                             | Age                              | 79       |      |
|                                 |                             | <input type="checkbox"/> Unknown |          |      |
| Study Information               |                             |                                  |          |      |
| Study ID                        | S00001                      |                                  | Year     | 2019 |
| Accession Number                | A00001                      |                                  | Month    | 07   |
| Study Description               |                             |                                  | Day      | 12   |
| Institution Name                | Hospital Name               |                                  | Hour     | 14   |
| Institution Address             | Hospital Address            |                                  | Minute   | 18   |
| Department Name                 | Department Name             |                                  | Second   | 28   |
| Referring Physician             |                             |                                  |          |      |
| Performing Physician's List     |                             |                                  |          |      |
| Operators' List                 |                             |                                  |          |      |
| Prefix                          | <input type="text"/>        |                                  |          |      |
| First Name                      | First                       |                                  |          |      |
| Middle Name                     | <input type="text"/>        |                                  |          |      |
| Last Name                       | Referring                   |                                  |          |      |
| Suffix                          | <input type="text"/>        |                                  |          |      |
| <div>OK</div> <div>Cancel</div> |                             |                                  |          |      |



## 5 Change items as desired.

Ensuring that the following value ranges are respected.

| No. | Group/Item          |                                     | Max Length | Permitted Characters/Numbers  | Max Total Characters |
|-----|---------------------|-------------------------------------|------------|---|----------------------|
| 1   | Patient Information | Patient ID                          | 64         | Alpha, Number, hyphen <->, period <.>, underscore <_>.  |                      |
|     |                     | Weight                              | 6          | Number, decimal point:<br>1 to 250 (kg),<br>2.2 to 550 (lb)   |                      |
|     |                     | Height                              | 6          | Number, decimal point:<br>30 to 242 (cm),<br>11.8 to 95.2 (inches)  |                      |
| 2   | Patient Name        | Prefix, First, Middle, Last, Suffix | 60         | <i>Combined length of all name fields cannot be greater than 60.</i><br>Alpha, Accents, hyphen, apostrophe <'>, underscore. | 60                   |
| 3   | Date of Birth       | Year                                | 4          | Number, no more than 107 years before current.  |                      |
|     |                     | Month, Day                          | 2          | Number. Month: 1-12, Day: 1-31.   |                      |
|     |                     | Age                                 | 3          | Number. 0-107.  |                      |
| 4   | Study Information   | Study ID                            | 16         | Alpha, Number, hyphen, period, underscore.  |                      |
|     |                     | Accession Number                    |            |   |                      |
|     |                     | Study Description                   | 64         | Alpha, Number, Accents, hyphen, period, apostrophe, underscore.   |                      |
|     |                     | Institution Name                    |            | Alpha, Number, Accents, Special.  |                      |
|     |                     | Department Name                     |            |   |                      |
|     |                     | Institution Address                 | 1024       | Alpha, Number, Accents, Special, backslash <\>.   |                      |
|     |                     | Year                                | 4          | Number. <i>Cannot be in future. No more than 107 years before current.</i>  |                      |
|     |                     | Month, Day                          | 2          | Number. 1-12, 1-31.   |                      |
|     |                     | Hour, Minute, Second                | 2          | Number 0-23, 0-59, 0-59.  |                      |
| 5   | Referring Physician | Prefix, First, Middle, Last, Suffix | 60         | <i>Combined length of all name fields cannot be greater than 60.</i><br>Alpha, Accents, hyphen, apostrophe.                 | 60                   |

| No. | Group/Item                  |                     | Max Length | Permitted Characters/Numbers        | Max Total Characters |
|-----|-----------------------------|---------------------|------------|-------------------------------------|----------------------|
| 6   | Performing Physician's List | First, Middle, Last | 60         | Alpha, Accents, hyphen, apostrophe. | 60                   |
| 7   | Operator's List             | First, Middle, Last | 60         | Alpha, Accents, hyphen, apostrophe. | 60                   |



The following accents are available:

Alpha: A-Z a-z

Number: 0-9

Accents: àâäçèéëêïïôöùûÛÀÂÇÉÊËËÎÔ <space>

Special: <backslash excluded>

- 6 To add a new performing physician, click [New] in the [Performing Physicians List], fill in the name fields, and click [Save]. To edit a performing physician, click the left or right-arrow button to display the previous/next name and then make the changes and click [Save]. To delete the displayed physician, click [Remove].
- 7 When finished making changes, click [OK] and then [Yes] to begin the modification process.  
Depending on the size of the study, it will take from several seconds to a minute or two to apply the modifications to every loop and image in the study. A progress bar grows toward the right during this process.
- 8 Once the modifications are complete, the Study List in the [Studies Management] window is refreshed, showing any modifications that you may have mad to the fields that appear in the Study List.

## 7.6 Anonymizing Study Information

A study can be copied and then have all patient-identity information replaced with generic values to protect a patient's identity. This is useful when images must be printed or transferred in a non-secure environment.

To create an anonymous copy or a study while leaving the original study untouched, follow this procedure.


- 1 Open the [Studies Management] window and click the study you wish to anonymize.
- 2 Click [Info]. The Stud Information Anonymization dialog box appears.
- 3 Click [Anonymize]. The Study Information Anonymization dialog box appears with Patient ID and Patient Name set to the generic value "Anonymous".

- 4 Make any additional modifications (respecting value ranges described in ["7.5 Modifying Study Information" P.7-12](#) and click [OK] and then [Yes] to begin the anonymization process.

Depending on the size of the study, it will take from several seconds to a minute or two to make a complete anonymized copy of the study including all loops and images. A progress bar grows toward the right during this process.

- 5 Once the anonymization is complete, the Study List in the [Studies Management] window is refreshed to show the just-added anonymous study.


By default, the new anonymous copy of the study retains the original accession number, however, its time and date are changed to the anonymous copy creation time.



## 7.7 Working with Image Annotations

Any area of interest in an image can be annotated with arrows and text superimposed on the image.

In this sample image, a stenosis has been highlighted with the annotation tool.




### NOTE

After saving an image, you can no longer annotate it.

To annotate an image, follow this procedure.


- 1 Pause a loop on the image you wish to annotate.

Click [Func]-[Abc Edit] on the side menu. The [Edit] tool bar appears floating over the image.



### 2 To add text, click [Text] button.

Move the pointer to where you wish to add text, click and begin typing. The text appears in the [Properties] dialog box and it is also superimposed over the image. Set color and style as desired.



### 3 Position the text by dragging it over the image. To add an arrow, click [Arrow] button, point to the area of interest in the image and click.

An arrow appears. Move the arrow around the image by dragging its head. Change the tail angle or length by dragging the tail. Set its color properties as desired.

### 4 If QCA or LV analysis feature is active while creating an annotation, the analysis information is superimposed over the image and the [Analysis] button on the [Edit] tool bar is automatically activated.

To hide this information, click the [Analysis] button so that it is not longer pushed in.




When image annotations are saved with the study as DICOM objects, all colors are converted to black and white. So once an annotated image has been saved with a study, closed, and then reopened, any colors will be gone. If desired, you can externally save the annotated image, including any color effects as described in ["7.10 Modifying Image Comments" P.7-33](#).

### 5 To add a copy of the Study Information displayed at the top of the Image Viewer window to the annotated image, click [Study] on the [Edit] tool bar.

A movable Study Information line appears over the image. Position it as desired.

- 6 When selecting other loop, the following dialog will display.



Click [No] to maintain the annotation.

- 7 Save your annotated image as described in the next section ["7.8 Saving Images" P.7-20](#).

### 7.7.1 Measurement Tools

ISO center value is used for calibration coefficient as a default. Once an image has been calibrated as described in ["8.4 Image Calibration Methods" P.8-11](#), three additional tools become available on the [Edit] tool bar. The [Caliper], [Line] and [Curve] tools can be used to make measurements in an image.

#### ■ Caliper

Select the [Caliper] tool and click on the image. The caliper appears. Grab the Caliper in the center and move it where you like. Size the Caliper by clicking the largest circle and then dragging inward to shrink it or outward to expand it.

#### ■ Line and Curve

Select the [Line] or [Curve] tool and click on the image. A point appears. Click elsewhere on the image to add a second point, and so on. The straight line (curve) connecting those points, is displayed. The dotted line is also displayed starting from the first point, and the length is displayed at the edge of the dotted line. If new points are added or deleted at the edge of the line, the length value is updated.

To stop entering points, click the [Pointer] button. To delete a point, select it at the edge of the dotted line that shows the straight line (curve) length and then click both mouse buttons simultaneously.

The edge point of the dotted line can be moved. When the point is moved, it is possible to change the position of where the straight line (curve) length is displayed. To move the point, drag it to a new position.

## 7.8 Saving Images




Save is only available for saving images from loops displayed on the Reference monitor.  
Save is not available for pre-existing still images.

The [Save] button on the side menu can be used to save still images (with any annotations) from loops when a study is open for acquisition or review.

To save an image, follow this procedure.

- 1 Pause on the image if it is a playing loop then click [Func]-[Save]. The following prompt appears.




- 2 Click [Yes] to perform the save. Any color in the image is converted to black and white before saving. Note that once you click [Yes], you cannot delete the new image object.

A new icon is added to the end of the Image Selector list for each image you save. The upper-left corner of the icon includes the loop number from which the image was saved, and in square brackets, the image frame number.

Also, when back in the [Studies Management] window, if you select this study and then click [Series], you will see your just-added object indicated by an increase in the number of objects in the [Annotated Image] series.



(Primary Archive mode only, for unprotected studies) When you save an image while in review, the entire study becomes protected because of the change to the study. See the Primary Archive section of  ["System Tab" P.17-17](#).



## 7.9 Print Images

### 7.9.1 Selecting Images for Printing

To select images for printing, including annotated and reference images, navigate to each image that you wish to print, pausing on it if it is in a playing loop, and click [Func]-[Add] on the side menu. Continue adding as many images as you like. When ready to see the print preview, click [Func]-[Print].



#### NOTE

The [Add] and [Print] side-menu buttons are grayed out if no printers are configured. Consult your system administrator for assistance.

If you click [Print], the print window appears like this.



You can then choose the desired DICOM printer in the Printer list. You can also configure Headers, Footers and Legends for both printer types. With the DICOM printer, you can print on special media such as photographic paper and film.

### 7.9.2 Printing Basics

---

After choosing [Funct]-[Add] to add images, click [Funct]-[Print]. A preview of your selected images is shown, arranged according to the Page setup values. Sufficient pages are added to accommodate the images. For example, if you add a total of 8 images and Columns=2 and Rows=2, (4 images per page), 2 pages will be required.

- 1 Navigate through the preview images with the [<<] and [>>] buttons in the lower-left corner.
- 2 Add images with the [Add] button.  
If desired, you can close the Print window to see the images better and then click [Funct]-[Print] again to reopen the window.
- 3 Delete images, by selecting a preview image (it is highlighted with a box) and then clicking [Delete].  
Click [Delete All] to delete every image from the print list.
- 4 In the Page setup box, set the number of Rows and Columns of images per page.  
For example, setting both Rows and Columns to 3 will give you 9 images per page.
- 5 Choose the desired printer in the Printer list.
- 6 Optionally configure printer features and set DICOM printer options as described in the two following sections.  
When ready, click [Print] to print the images.
- 7 If desired, you can click [Close] without printing and print later.  
Your images are preserved in the Print list until you delete them or close the study.



Printed images can be of diagnostic quality only when printed on an approved medical-grade printer that has been calibrated.

### 7.9.3 Print Features

Print Features such as Headers and Footers can be configured.



To print in color (if possible), all color information added to the image will be printed in color. However, only black and white images will be saved, but not color images.

#### Horizontal/Vertical

In the Page setup box, you can choose [Horizontal] to fill the page left to right and then down, or [Vertical] to fill it top to bottom and then to the right. For example, if Columns and Rows are set to 2, for 4 images per page, selecting [Horizontal] will cause the first image to be placed in the upper-left area, the second to be placed in the upper-right area, the third to be placed in the lower-left area, and so on. In contrast, [Vertical] will cause the second image to be placed in the lower left and the third image to be placed in the upper-right.

#### Modifying the Header and Footer

By default, a header and footer is included on every page. You can use the defaults, adjust them, or completely replace them.

To modify the Header or Footer, follow this procedure.

- 1 Click the [Header] or [Footer] button.

The Header or Footer edit box appears.

**Header**

|                           |                   |
|---------------------------|-------------------|
| Institution Name          | :%InstitutionName |
| Patient Name              | :%PatientName     |
| Performing Physician Name | :%PerPhysName     |
| Accession Number          | :%AccessionNumber |
| Patient ID                | :%PatientID       |
| Study Date                | :%StudyDate       |

Font: B I U Alignment: ☒ Left ☐ Right ☐ Center Institution Name Add Save Close

**Footer**

Page %P of %N

Font: B I U Alignment: ☐ Left ☒ Center ☐ Right Institution Name Add Save Close


- 2 Enter any text you like. Press [Enter] key to start a new line.  
Press the keyboard [Tab] key to align text in columns. To delete text, select it (drag over text with mouse button pressed so that selection becomes highlighted) and press the keyboard [Delete] key.


- 3 Variables and macros for common information such as Patient Name can be included by selecting them from the drop-down list above the [Add] button and then clicking [Add].

During preview and printing the variables are replaced with the actual values from the study. You can add, delete, and rearrange macros and variables or your own text in any way desired. The variables provide information specific to the image added to the print list.

The following variables are supported:

| Variable  | Notes  |
|---|--|
| General: These variables provide general information about the study.                               |  |
| %InstitutionName  |  |
| %PerPhysName  |  |
| %PatientName  |  |
| %StudyDate  |  |
| %PatientID  |  |
| %StudyID  |  |
| %Patient Sex  |  |
| %Patient Birth Date   |  |
| %AccessionNumber  |  |
| %Page Number  |  |
| %Number of pages  |  |
| Parameters: These variables replicate any parameter information associated with the selected image. |  |
| %Oblique Angle  | Oblique angle of the detector in degrees in which angles toward the patient's left are displayed with LAO, and angles toward the patient's right are displayed with RAO. |
| %Sagittal Angle   | Sagittal angle of the detector in degrees in which angles toward the patient's head are displayed with CRA, and angles toward the patient's foot are displayed with CAU. |
| %sid  | X-ray source to detector distance in centimeters.  |
| %fov  | FPD field of view in inches.   |
| %kV   | X-ray tube voltage in kilo-volts.  |
| %mA   | X-ray tube current in milliamperes.  |
| %mAs  | Function of X-ray tube current and time expressed in milliamperes second.  |
| %Time   | X-ray pulse (exposure) duration for this image, in milliseconds.   |
| %Image Comment  | The selected comment (if any).   |
| %Dose Of The Loop   | Total dosage of the loop.  |

| Variable   | Notes   |
|--|---|
| %Image Number  | Number of the image.  |
| %Acquisition Time  | The time which the image was taken.   |
| Frame Information: These variables provide frame information.  |   |
| %Frame   | The image from which this image was taken.  |
| %n_frame   | The total number of frames in the loop.   |
| %Frame Time  | The number of seconds since start of loop at which the image was acquired.  |
| %Injection Time  | The period of time when turn ON the injector with overlay by frame in milliseconds.   |
| %Inverted  | The image is inverted.  |
| Quantitative Coronary Analysis (QCA):  <a href="#">"8.2 Quantitative Coronary Analysis" P.8-2</a> |   |
| %QCA Reference   | <p>These variables are used inside the %qca macro to provide QCA result values for images with active QCA analysis at the time they are added to the print list.</p> <p>The default <b>%qca</b> macro definition is as follows:</p> <pre>%qca( QCA ref %QCA Reference Stenosis %QCA RAO Stenosis mm mm : %QCA RAO Steno- sis % % )</pre> <p>Enter this entire definition in one of the legend boxes. This causes QCA results to be displayed like this:</p> <div style="background-color: black; color: white; padding: 5px; text-align: center;"> <b>QCA ref 3.96</b><br/> <b>Stenosis 1.43 mm : 63.84 %</b> </div> <p>The variables in this macro only produce results for images on which QCA results had been calculated at the time the image was added to the print list. If this is not the case, no QCA result values are printed.</p> <p>If desired, you can omit some of the individual variables or re-sequence them.</p> <p>Although not used in the default %qca macro definition, variable %Calib Type provides calibration type information like this, according to selected calibration type:</p> <p>Distance: 20.0 mm<br/> Grid: 2x2<br/> Catheter 10F 3.33 mm</p> |
| %QCA RAO Stenosis mm   |   |
| %QCA RAO Stenosis  |   |

| Variable   | Notes  |
|--|--|
| Left Ventricular Analysis (LVA):  "8.3 Left Ventricular Analysis" P.8-7 |  |
| %LV RAO EDV ml   | <p>These variables are used inside the default <b>%lvrao</b> macro to provide LV RAO values for images with active LV analysis at the time they are added to the print list.</p> <p>The default <b>%lvrao</b> macro definition is as follows:</p> <pre>%laver( RAO 30° Defame) : %LV RAO EDV ml Defame) : %LV RAO EDV m2 ESV(m2) : %LV RAO ESV ml ESV(m2) : %LV RAO ESV m2 EF : %LV RAO EF</pre> <p>Enter this entire definition in one of the legend boxes. This causes LV RAO results to be displayed like this.</p> <div data-bbox="912 777 1134 1075" data-label="Text"> <pre>RAO 30° EDV(ml) : 65.19 EDV(m²): 27.34 ESV(ml) : 14.15 ESV(m²): 6.16 EF : 79.40</pre> </div> <p>If desired, you can omit some of the individual variables or re-sequence them.</p> <p>The patient's height and weight must be entered on the LV tool bar before proper result can be calculated.</p> |
| %LV RAO EDV M <sup>2</sup>   |  |
| %LV RAO ESV ml   |  |
| %LV RAO ESV M <sup>2</sup>   |  |
| %LV RAO EF   |  |

| Variable                       | Notes  |
|--------------------------------|--|
| %LV RAO LAO EDV ml             | These variables are used inside the default <b>%lvraolao</b> macro to provide LV RAO values for images with active LVA analysis at the time they are added to the print list.  |
| %LV RAO LAO EDV m <sup>2</sup> |  |
| %LV RAO LAO ESV ml             |  |
| %LV RAO LAO ESV m <sup>2</sup> |  |
| %LV RAO LAO EF                 | <p>The default <b>%lvraolao</b> macro definition is as follows:</p> <pre>%lvraolao( RAO 30°+ LAO 60° EDV(ml) : %LV AO LAO EDV ml EDV(m2): %LV RAO LAO EDV m2 ESV(ml) : %LV RAO LAO ESV ml ESV(m2) : %LV RAO LAO ESV m2 EF : %LV RAO EF</pre> <p>This causes LV RAO LAO results to be displayed like this:</p> <div style="background-color: black; color: white; padding: 10px; text-align: center;"> <p>RAO 30° + LAO 60°</p> <p>EDV(ml) : 45.06</p> <p>EDV(m<sup>2</sup>): 18.99</p> <p>ESV(ml) : 5.87</p> <p>ESV(m<sup>2</sup>): 2.73</p> <p>EF : 79.40</p> </div> <p>If desired, you can omit some of the individual variables or re-sequence them.</p> <p>The patient's height and weight must be entered on the LV tool bar before proper results can be calculated.</p> |



Special variables, such as those containing QCA or LV analysis results (e.g., %QCA RAO Stenosis %), will only have a value for images that were added to the print list when the QCA or LVA results were displayed.

- 4 To format, select the desired text and variables, and then click one or more of the buttons: [B] for bold, [I] for italic, [U] for underlined, or [Font] for font.  
For Font, choose Font, font style, and Size from the dialog box that appears.
- 5 Change text alignment by selecting one or more lines of text and clicking an Alignment button: [Left], [Right], or [Center].
- 6 When finished, click [Save]. Your changes are immediately reflected in the Preview.



## ■ Modifying the Image Legend

Any Legend text you like can be superimposed over the four corners of the image. To add Legend text, click the [Legend] button, enter your text and variables, and then apply formatting like done for Header in any of the four boxes.

**Legend**

**Top Left**

%qca(  
QCA ref %QCA Reference  
Stenosis %QCA RAO Stenosis mm mm : %QCA RAO Stenosis % %  
)  
%lvrao(  
RAO 30°  
EDV(ml) : %LV RAO EDV ml  
EDV(m²) : %LV RAO EDV m2  
ESV(ml) : %LV RAO ESV ml

**Top Right**

**Bottom Left**

%Oblique Angle : %Sagittal Angle  
SID %sid  
%fov inch  
%kVpV %mA %Time

**Bottom Right**

%Inverted  
%Frame/%n\_frame  
%Frame Time

**Font** **B** **I** **U**

**Alignment**  
☒ Left ☐ Right  
☐ Center

%Oblique Angle

**Add** **Save** **Close**

Click [Save] when done. Add a new image and check the result in the Preview. Until modified, this legend will be included on all images subsequently added to the Print list.

## ■ Selecting Page Range

If you wish to print fewer than all image, select the desired page(s) in the Page range box.

## Setting Printer LUT

If configured at your site, you may be permitted to select alternative printer LUT profiles to control image brightness and contrast. If the current printer LUT profile (name displayed near LUT button) is not appropriate, follow this procedure to set a different one.

### 1 Click [LUT] button.

A small profile selection box appears at the upper-right corner of the [Print] window.



### 2 Choose the desired profile from the drop-down list.

The chosen profile name appears near [LUT] button. If desired, click [Close] to close the drop-down list.


## Setting Printer

Select from [Printer] on the [Print] window and click [Properties]. The following [Printing Configuration] window appears.


### 7.9.4 Setting DICOM Printer

There are several DICOM printer properties that are configured as follows.

Select [Funct]-[Print] from the side menu and click [Properties] for [Printer].



The following [Printing Configuration] window is displayed.




| No. | Item                            | Description   |
|-----|---------------------------------|---|
| 1   | Orientation                     | Choose [LANDSCAPE] for a page that is wider than it is long, or [PORTRAIT] for the opposite.  |
| 2   | Media Type                      | Choose media such as paper, photographic paper, or film, according to what is supported and available. Black text on a white background is used for paper, and white (clear) text on a black background is used for film. |
| 3   | Magnification                   | Choose the scaling algorithm used to fit images into the designated page region. For example, BICUBIC will produce better quality than CUBIC, however it will take more time.   |
| 4   | Polarity                        | Only for the image itself, choose [NORMAL] for normal printing or [REVERSE] for reverse (white inverted to black, black inverted to white) printing.  |
| 5   | Smoothing                       | Select smoothing level of the image from a pull-down list. Enabled only for some [Magnification] values (e.g., [CUBIC]). Select a value supported by the selected printer.  |
| 6   | Media Size                      | Choose media size.  |
| 7   | Max Horizontal Size             | Select the document maximum horizontal size (in inches).  |
| 8   | Max Vertical Size               | Select the document maximum vertical size (in inches).  |
| 9   | Get Printer Status / Get Status | Click to display DICOM printer [Status].  |


## 7.10 Modifying Image Comments

Enable to modify comments on the image after study. Follow the procedure below to modify the comments.

- 1 Display an image which comment to be modified.
- 2 Click [Comment] on the side menu.
- 3 Select [New Comment] to set a new image.
- 4 After entering a comment, click [Modify] button.



- 5 To select from registered image comments list, select [Comments List].



- 6 After selecting an image comment, click [Modify] button.

This page is intentionally left blank.

# Chapter 8

## Performing Angiographic Analysis

This chapter describes the Quantitative Coronary Analysis (QCA) and Left Ventricular (LV) analysis features and their related calibration procedures. This chapter provides detailed information on the analysis process as follows.

### Description

|     |                                      |      |
|-----|--------------------------------------|------|
| 8.1 | Introduction .....                   | 8-2  |
| 8.2 | Quantitative Coronary Analysis ..... | 8-2  |
| 8.3 | Left Ventricular Analysis .....      | 8-7  |
| 8.4 | Image Calibration Methods .....      | 8-11 |

## 8.1 Introduction

Quantitative Coronary Analysis (QCA) and Left Ventricular (LV) features are available to help assess coronary stenosis and ventricular volumes.

## 8.2 Quantitative Coronary Analysis

The QCA feature enables you to accurately measure coronary artery stenosis. In automatic mode, you just click a few points in the general area of the stenosis and the QCA feature does the rest, drawing a contour and calculating the degree of stenosis.




### WARNING

Measuring a vessel diameter from angiography causes several errors. Therefore, do not diagnose or treat based solely upon this analysis result. Use the analysis result as a reference value.


To measure coronary artery stenosis, follow this procedure.

- 1 Display a study containing both calibration reference and stenosis images. Click [Funct]-[QCA].

The image is automatically zoomed and the [QCA] tool bar is displayed. Calibration and analysis information is superimposed over the image, if configured (see  "Display Tab" P.17-37).






- 2 ISO center value is registered for calibration coefficient as a default value. For more detailed analysis result, perform calibration as described in  "8.4 Image Calibration Methods" P.8-11.  
This determines the mm to pixel ratio of the image. When calibrated, display an image with a good view of the stenosis.
- 3 If needed, pan (right-click drag) the image so that the area of interest is centered in the Image Viewer.
- 4 You can also pan the image by dragging the pan box at the bottom of the tool bar with either the mouse (left-click drag) or, if configured, the joystick (hold down button)  
Drag the pan box in the direction you wish to pan the image.
- 5 You can zoom in and out the displayed image.  
Right-click on the following button on the side menu to zoom in and left-click to zoom out.




- 6 Click the [Automatic]-[Stenosis] button. Then click three points along the stenosed artery.

With the first point on one side of the stenosis, the second point near the stenosis, and the third on the other side of the stenosis. Then click a fourth point on the part of the artery you judge to be the coronary reference diameter.




- 7 A contour is drawn along the artery.

With the reference point marked by a white dot with a line through it extending to the artery edges. the narrowest part of the artery, the suspected stenosis, is marked with a white line across the artery. The degree of stenosis is calculated based on the ratio between the width of the artery at the reference point versus the width of the artery at the most stenosed point. Drag the reference line along the artery, or adjust other points as desired. The percent stenosis is recalculated instantly.




- 8 To adjust an auto-generated contour line to fit the artery edge more closely, click [Edit] on the [QCA] tool bar.

The image is zoomed in so you can see the artery more closely. Right-click drag the image to center the area of interest. Click [Redraw] on the Edit tool bar.




- 9 On the image, click where you want to start adding points and then click to add other new points along the artery edge where the automatic contour is not close enough. Drag the new points as needed to follow the artery edge. The first and last new points automatically connect to the automatic contour. Click [Apply] to change the contour and recalculate stenosis.

When adjusting the contour, press the keyboard space bar to temporarily hide the automatic contour so that you can clearly see the artery edge.




- 10 The contour is modified, the stenosis indicator line is moved to the new narrowest point.


Stenosis is recalculated and displayed.




- 11 To trim off the ends of the contour lines, click [Edit] on the [QCA] tool bar and then [Cut] on the [Edit] tool bar.



- 12** Click points on either side of the ends of the contour lines like this so that the connecting line bisects the contour lines at the desired trim points.




- 13** Click [Apply] to trim the contour lines where indicated.




- 14** To measure stenosis manually, click [Manual]-[Stenosis] button. Click points on opposite sides of the narrowest part of the stenosed artery.

The stenosis is measured and indicated on screen.



- 15** Click [Manual]-[Ref.] button. Click points on opposite sides of the part of the artery that you consider to be un-stenosed.

The reference distance is measured and indicated on screen. The percent stenosis is also calculated and displayed.



- 16** Select [Line] or [Curve] and click on the image to show points. Lines or curves connecting the points are displayed.

A dotted line which begins from the first point is displayed at the same time, and length is displayed at the end of the dotted line. The length of line and curve displayed at the end of the dotted line is updated when adding or deleting the points.

- 17 To delete the points, select the end of dotted line showing the length of straight line or curve, and click both right and left mouse buttons at the same time.
- 18 Click [Delete] to delete currently selected object.
- 19 Click [Del. All] to delete all the objects except QCA calibration object.

## 8.3 Left Ventricular Analysis

The left ventricular analysis feature enables you to calculate left ventricular volume and analyze the regional and global wall motion with the Modified Centerline method. These values are determined starting with a right oblique anterior projection at 30 degrees, a left oblique anterior projection at 60 degrees, and their spatial combination.



### CAUTION


For optimal Left Ventricular analysis, both the heart and the calibration device should be positioned at the center of the image for RAO 30 degrees and LAO 60 degrees. In other words, they should be positioned at ISO center.


To perform left ventricular analysis, follow this procedure.

- 1 Display a study containing both a calibration reference image and other images showing the left ventricle in its various states.

- 2 Click [Funct]-[LV].

- 3 The [LV] tool bar is displayed.

Calibration and analysis information is superimposed over the image, if configured. See  ["Display Tab" P.17-37](#).



- 4 If needed, pan (right-click drag) the image so that the area of interest is centered in the Image Viewer window.

You can also pan the image by dragging the pan box at the bottom of the tool bar with either the mouse (left-click drag) or the joystick (hold down button). Drag the pan box in the direction you wish to pan the image.

- 5 ISO center value is registered for calibration coefficient as a default value. For more detailed analysis result, perform calibration as described in ["8.4 Image Calibration Methods" P.8-11](#). This determines the mm to pixel ratio of the image.



If the image has not been calibrated, only the results of the ejection fraction and segmental LV function will be available.

- 6 In the Contour group of the [LV] tool bar, click the projection icon that corresponds to the image (the End of Diastole image (ED), or the End-of Systole image (ES) of RAO or LAO projection).



RAO ED



RAO ES




LAO ED




LAO ES

- 7 Working with only RAO ED and RAO ES projections, two contours can be traced, and if the LAO projection is also used, four contours may be traced.  
The selected projection is indicated by the highlighted projection icon on the LV tool bar.
- 8 Using the mouse, outline the ventricle's contour by clicking multiple points along the ventricular wall. Start contouring from the base of the aorta.  
It is sufficient to click a few points along the contour by clicking the left mouse button and then the outline curve is automatically created.
- 9 Move a point by clicking and dragging it to a new position.
- 10 Delete a point by clicking it with the left and right mouse buttons simultaneously.
- 11 To erase the entire curve, click [Clear]. Alternatively, to cancel recent changes and restore previous values, click [Cancel].
- 12 Once the contour is defined up to the lower-right aortal base, validate the contour by clicking the icon of the next projection to be drawn.  
Once validated, an opaque square surrounds the projection icon to indicate that the corresponding contour has been entered and validated.

- 13 Select the image for ES and define its contour.  
Here is a sample with the ED and ES contours completed.



- 14 Click [Compute] icon.  
The results of the ventricular analysis are displayed like this.



Abbreviations:

BA= Basal Anterior, AA= Apical Anterior, API= Apical, AI= Apical Inferior, BI= Basal Inferior, BL= Basal Lateral, AL= Apical Lateral, AS= Apical Septal, BS= Basal Septal



## CAUTION

This analysis is calculated with an approximation formula called “Area Length,” so the calculation result includes approximation errors.



### CAUTION

Analyze both LAO and RAO for high precision analysis. Also, Perform either of calibration once.



### CAUTION

Calibrate on either the LAO or RAO image. It is not necessary to calibrate on both. Be sure to adjust the target region to the ISO center and perform LAO and RAO radiography without changing SID. Also, perform either of calibration once.



## 8.4 Image Calibration Methods

Calibration options available for the QCA and LV features are as follows.



### CAUTION

If you choose a different loop after calibrating, a message may be displayed asking if you wish to keep the existing calibration. Keep the existing calibration **ONLY** if you know that the acquisition configuration and C-arm parameters are the same for the different loop.

### 8.4.1 Using ISO Center Value

ISO center value is used as a default, however, perform calibration by manual if needed.

### 8.4.2 Using a Catheter (QCA)

To calibrate an image using a catheter, follow this procedure.

- 1 Display an image that contains a catheter of known diameter.
- 2 On the [QCA] tool bar, click the arrow button below [Calibrate] and choose [Catheter]. Then click the [Catheter] button.
- 3 Choose one of the predefined catheter diameter values from the drop-down list or enter the exact value (mm) in the box provided.  
Enable to enter the value up to 50 mm.
- 4 Click [Catheter].
- 5 Click 3 points along the middle of the catheter.
- 6 The calibration factor in mm pixel is calculated and displayed in the [Calibration Factor] box.  
To adjust the outline and improve the mm per pixel accuracy, click any of the three points and drag it while keeping the mouse button pressed.



### CAUTION

If the target vessel for measurement and geometric expansion of catheter for calibration is different, then the expansion difference causes errors.

### 8.4.3 Using Reference Distance (QCA)

---



### CAUTION

This method is only valid if the object of known dimensions lies in the same plane as the structure to be measured in the image.

To calibrate an image using a reference distance, follow this procedure.

- 1 Display an image that contains an object of known dimensions.
- 2 On the [QCA] tool bar, click the arrow button below [Calibrate] and choose [Distance]. Then click [Distance].
- 3 Enter the known reference distance (mm) in the box provided.
- 4 Click [Distance] button.
- 5 Click two points on the image to mark the exact starting and ending points of the object of known dimensions. A line segment is drawn between the two points.
- 6 The calibration factor in mm per pixel is calculated and displayed in the [Calibration Factor] box.
- 7 If desired, adjust the line segment as follows.
  - Move the segment on the image by clicking the first point of the segment and dragging it while keeping the mouse button pressed.
  - Stretch or reduce the segment by clicking the last point of the segment and dragging it while keeping the mouse button pressed.

#### 8.4.4 Using a Grid (QCA and LV)



### CAUTION

This method is only valid if the grid of known dimensions lies in the same plane as the structure to be measured in the image.

To calibrate an image using a grid, follow this procedure.

- 1 Display an image that contains a grid of known dimensions.
- 2 On the [QCA] or [LV] tool bar, click the arrow button below [Calibrate] and choose [Grid]. Then click [Grid].
- 3 Choose a grid from the drop-down list, from 1 mm square up to 7x7 mm.
- 4 Click [Grid].
- 5 Click four points on the image to mark the initial position and size of the grid. A grid is drawn over the image. Now position the grid over a corresponding number of grid cells in the image.  
  
For example, if you chose a 2x2 grid, position the on-screen grid so that it is superimposed over a 2x2 grouping of grid cells in the image.
- 6 The calibration factor in mm per pixel is calculated and displayed in the [Calibration Factor] box.

#### 8.4.5 Using a Reference Ball (LV)

To calibrate an image using a reference ball, follow this procedure.

- 1 Display an image that contains a reference ball of known diameter.
- 2 On the [LV] tool bar, click the arrow button below [Calibrate] and choose [Ball]. Then click [Ball] again.
- 3 Enter the diameter of the ball (mm) in the box provided.
- 4 Click [Ball].
- 5 Click 3 points on the ball circumference. The ball is outlined and the calibration factor in mm per pixel is calculated and displayed in the [Calibration Factor] box.  
  
If desired, adjust the ball outline by clicking and dragging any of the 3 points.

This page is intentionally left blank.

# Chapter 9

## DSA Image Adjustment Tools

This chapter describes how to adjust Digital Subtraction Angiography (DSA) images with the tools provided.

### Description

|     |  |      |
|-----|--|------|
| 9.1 | Introduction . . . . .                     | 9-2  |
| 9.2 | Adjusting Subtraction Parameters . . . . . | 9-3  |
| 9.3 | Reregistration. . . . .                    | 9-5  |
| 9.4 | Peak Hold . . . . .                        | 9-6  |
| 9.5 | Contour Enhancement. . . . .               | 9-8  |
| 9.6 | Mask Region . . . . .                      | 9-11 |
| 9.7 | Adjusting Shutter . . . . .                | 9-13 |

## 9.1 Introduction


Images captured in DSA mode have constant structures of no diagnostic interest virtually removed, enabling enhanced blood-vessel contrast.

Tools are provided to adjust DSA images to ensure the best-possible image quality.


## 9.2 Adjusting Subtraction Parameters


To modify the subtraction parameters of a loop.

- 1 Display the DSA loop.
- 2 Click [Funct]-[DSA] on the side menu.
- 3 The [DSA Tools] bar is displayed with the frame numbers of the mask frames shown in the [Mask #] box. The default number of consecutive frames used to create the mask (defined in Rad program definition) is shown in the [Mask #] selection box.



- 4 To change from the automatically-selected mask frames, pause on the image frame to be used as the first frame for the new mask and click [Set Mask] button on the [DSA Tools] bar. This frame and the appropriate number (defined in Rad program definition) of additional consecutive frames are used as the mask.
- 5 To change the number of additional consecutive frames used to create the mask, select a different value from the [Mask #] drop-down list of the [DSA Tools] bar.
- 6 To change the Mask Weight, use [Mask Weight] button to change the value (from 0 to 100%).



Increment for the Up/Down button can be changed in the range of 1 to 25.  ["Display Tab" P.17-37.](#)

After changing the increment, close the [DSA Tools] bar, and then open it again.

- 7 If necessary, display the image in subtracted mode and use the side menu buttons to adjust the gray scale.

- 8 For DSA images, subtraction mode is the default. To review images in non-subtraction mode, click [Subtract] button on the [DSA Tools] bar to deselect it (no longer appears pushed in).  
Click [Subtract] again to return subtracted mode.
- 9 To change blood vessel contrast in non-subtraction mode, change the value of [Vessel Power] (from 0 to 1000).



## 9.3 Reregistration

If necessary, use the [Pixel Registration] tool to reposition the image mask.



- 1 Click [Shift] on the [DSA Tools] bar.
- 2 Move the image along any of the four-axis by clicking the appropriate arrow button on the [Pixel Registration] tool bar.
- 3 If it is desired to move less than one pixel per click, select a different interval increment from the [Interval] drop-down list on the [Pixel Registration] tool bar.
- 4 To reposition the mask to the initial position (0,0), click [Reset] on the [Pixel Registration] tool bar.
- 5 Acquired with Flex-APS (option) and position adjusted loop is checked on the check box of [Flex-APS]. Uncheck this box to stop position adjustment.

## 9.4 Peak Hold

Create peak hold image.




**NOTE**


Cannot use Mask Weight in combination with Peak Hold.

- 1 Click [Peak Hold] button on the [DSA Tools] bar and the peak hold will automatically apply to displaying image.

The peak hold setting is that peak mode is black and the frame range is currently selected frame to the last frame.



- 2 In case of that adjustment of peak hold image is necessary, use the peak hold dialog and set as follows.



- 1 Select either black or white peak.
- 2 Frame range applying peak hold. Display Start frame and click [Set Frame] button. And then display [End frame] and click [Set Frame] button.
- 3 Click [Compute] to recreate the peak hold image.

Compute

#### Using Mask Image

Enable to use created Peak Hold Image as a mask image.

To use as a mask image, click [Fluoro]-[DSA MAP].

- 4 Click [Save] to save the created peak hold image as a still image.

Save



Still image cannot register as a mask image. To superimpose the image on fluoroscopy image, use the image just after creating the peak hold image.


## 9.5 Contour Enhancement

Use contour function to create an contour enhanced blood-vessel image from DSA image.  
Superimpose and display the created image on the fluoroscopy image.



Contour Enhancement is a function to know the shape of blood-vessel and branch position for large blood-vessel diameter such as aorta. Therefore, when apply this function to small blood-vessel, the blood-vessel may not be rendered, or rendered blood-vessel diameter may be too wide or too narrow.


- 1 Display the target image on REF monitor.
- 2 Click [Contour] button on the [DSA Tools] bar.



- 3 The image on the monitor will be edge enhanced with default parameters.



- 4 In case of that adjustment of vessel contour is necessary, set parameters of [Region Rendering Quality], [Smoothing] and [Thickness] on the Contour Enhancement dialog.



Description of each parameter is shown in the following table.

|   |                          |   |
|---|--------------------------|---|
| 1 | Region Rendering Quality | Click [+1] or [+10] button to render wider region of blood-vessel for displaying more less density of blood-vessel contour. And click [-1] or [-10] button to render narrower region of blood-vessel, if rendered region is too wide or too much noise. |
| 2 | Smoothing                | Enable to specify the smoothness of displaying blood-vessel. Click [Smooth] button to smooth the blood-vessel image and [Sharpen] to reflect accurate shape of the blood-vessel.  |
| 3 | Thickness                | Enable to change contour thickness.   |


Enable to set default parameters for each physician. "Devices Tab" P.17-56


- 5 Click [Save] button to save the image as a still image.



#### NOTE

Still image cannot register as a mask image. To superimpose the image on fluoroscopy image, use the image just after applying Contour processing.


- 6 After creating the contour enhancement image, click [MAP]-[DSA] button on the side menu. Enable to register the image as a Mask image.  ["DSA-MAP" P.4-32](#)



## 9.6 Mask Region

To superimpose only a part of blood-vessel image of peak hold, contour enhancement and DSA images, enable to set ROI of SIMAP.

- 1 Display an image such as peak hold, edge enhancement and DSA images to be registered as SIMAP.
- 2 Click [Region] button on the [DSA Tools] bar.




- 3 ROI is displayed on the image.



- 4 Set ROI where you want to use as MAP image on the displayed image. Drag a mouse on ROI to move and drag the edge of ROI to modify the size. Enable to move and modify the size with IVR NEO/SMART Touch.
- 5 Click [MAP]-[DSA] button on the side menu. ROI area is superimposed and displayed on subsequent fluoroscopy image.

 ["DSA-MAP" P.4-32](#)





## 9.7 Adjusting Shutter

Shutter completely block out image edges. Use shutters to hide distracting portions of an image, so that you can focus on areas of interest.

Combined circular and rectangular shutters are always adjustable. Shutter positioning and sizing information is retained for each group.

to move the shutter edges, position the mouse pointer at the edge of the image and find the shutters which are displayed as a white rectangle and circle. Use the left mouse button to resize the shutters by dragging a shutter edge inwards to make the visible image area smaller or outwards to make it larger.

The four sides of the rectangular shutters can be individually dragged to their desired position. The image shows through the intersecting portion inner area of the combined circular and rectangular shutters.

This page is intentionally left blank.

# *Chapter 10*

---

## Difference with Bi-plane Option

This chapter describes differences when using the Bi-plane option to provide simultaneous frontal and lateral imaging.

### ■ Description

|      |                                  |       |
|------|----------------------------------|-------|
| 10.1 | Introduction .....               | 10-2  |
| 10.2 | Bi-plane Image Acquisition ..... | 10-5  |
| 10.3 | Bi-plane Image Preview .....     | 10-14 |
| 10.4 | Shutdown .....                   | 10-15 |
| 10.5 | Administration .....             | 10-16 |

# 10.1 Introduction

## 10.1.1 Standard Components

The standard components of a bi-plane system are as follows.

| No. | Name  | Installation Location | Quantity |
|-----|---|-----------------------|----------|
| 1   | Digital reference cabinet:<br>With built-in Reference computer                                    | Control room          | 1 set    |
| 2   | Digital acquisition cabinet:<br>With built-in acquisition computer, CP2 (Panel Control Unit), UPS | Equipment room        | 2 sets   |
| 3   | FPD cooling unit  | Equipment room        | 2 sets   |
| 4   | Monochrome LCD monitor  | Control room          | 4 sets   |
|     |   | Examination room      | 4 sets   |
| 5   | Keyboard and mouse  | Control room          | 1 set    |
| 6   | Bedside console IVR NEO   | Examination room      | 1 set    |
| 7   | FPD (Flat Panel Detector)   | Examination room      | 2 sets   |
| 8   | FPD high and low-voltage power unit   | Examination room      | 2 sets   |
| 9   | IVR Shuttle (Option)  | Control room          | 1 set    |


## 10.1.2 Bi-plane User Interface

With the bi-plane option, a second C-arm, X-ray generator and acquisition computer are added for the Latest plane so that simultaneous Frontal/Lateral imaging can be performed. A second Reference monitor, the Reference-Lateral monitor, is included for synchronized Frontal and display Lateral image.

Images and loops acquired on a Bi-plane system include a plane indicator in the icon labels shown in the Image Selector (A for plane A (Frontal) and B for plane B (Lateral)). The Plane A/B designation is also included in the DICOM tags. An active-plane indicator icon is displayed in the top left corner of the Reference monitors.


For example, the Reference-Frontal monitor showing the Frontal images and the Reference-Lateral monitor showing the Lateral images.

- Reference-Frontal Monitor



For detailed descriptions of these window elements, see [☞ "4.2 The Main User Interface Window" P.4-2.](#)

- Reference-Lateral Monitor



Loops 1A and 1B were acquired at the same time, loops 2A and 2B were acquired at the same time, and so on. Additionally, the image frames are acquired in synchronization, meaning that 1A frame 1 corresponds to 1B frame 1, 2 to 2, 3 to 3, and so on.



When using a Bi-plane system as a Single-plane system, move the Lateral C-arm to the park position and use the Frontal arm only.


### 10.1.3 Reference Monitors

Only the Reference-Frontal monitor displays the full GUI with side menu and bottom bar. Both Reference monitors display their respective (Frontal or Lateral) loops, image, Image Selectors, identification and parameters. Frontal/Lateral Loop playback is always synchronized.


Mono-plane loops (loops acquired on a bi-plane system but only on one plane) are displayed on their respective Reference monitor with the other Reference monitor paused, showing an "X" mark in the upper-left check box. The Frontal Image Selector shows only Frontal image and loop icons. The Lateral Image Selector shows only Lateral image and loop icons.

# 10.2 Bi-plane Image Acquisition

## 10.2.1 Entering Study Information


Study information is entered in the same way as for non-bi-plane systems. Follow the procedure in  "4.3 Entering Study Information" P.4-10.

## 10.2.2 Acquiring Loops when Bi-plane Active

When Bi-plane is active, loops are acquired in similar manner to when it is not active (see  "4.4 Making Acquisition" P.4-23).

- Fluoro acquisition is performed on the Frontal or Lateral planes according to the Fluoro switch pressed. Press either the Frontal Fluoro switch or the Lateral Fluoro switch.
- When pressing both Frontal and Lateral fluoro switches at the same time at 30 pps, you cannot perform Bi-plane fluoroscopy (irradiating the Frontal and Lateral planes at the same time).
- Rad acquisition is performed on the Frontal, Lateral or simultaneously on both planes, according to the program selected. Choose Rad programs as follows.

1 Click [Program: Rad] or press the button on IVR NEO/IVR Shuttle/SMART Touch.




Side Menu



IVR NEO





IVR Shuttle




SMART Touch

2 Select Rad program.



- 3 To change the plane, click [Program: Change Acquisition Plane] and select the desired plane. Currently selected DUP can be changed to F, L or Bi-plane.



On DUP which acquisition plane is restricted, disable to change the plane with [Change Acquisition Plane] button.

- Images from one or both planes are then transmitted to the Reference system. Bi-plane images are displayed synchronized on the two monitors.
- For DSA-MAP on Bi-plane, enable to register the following images as a blood-vessel image depending on the monitor selecting [DSA-MAP] button on the side menu.
  - 1 When selecting on the Frontal Acquisition monitor, register the displayed image on this monitor as a blood-vessel image.
  - 2 When selecting on the Acquisition monitor, register the displayed image on this monitor as a blood-vessel image.
  - 3 When selecting on the Reference monitor and if it is in RefSync mode, register the image on both Frontal/Lateral Reference monitors as blood-vessel images at the same time on each Acquisition monitor. And if it is not in RefSync mode, register the selected plane image as a blood-vessel image.
- When selecting DSA-MAP on the Reference monitor, if it is in RefSync mode, register the image on both Frontal/Lateral Reference monitor as blood-vessel images at the same time. And if it is not in RefSync mode, register the selected plane image as a blood-vessel image.



### 10.2.3 Bi-plane Loop

#### Reference Monitor

- **Ref Synchronized Mode**

In Ref Synchronized mode (RefSync side menu button pushed in), image playback (loop and frame) direction and speed (fps), and zoom are always synchronized on the two Reference monitors.



Other functionality is controlled separately on the selected plane (Reference-Frontal or Reference Lateral.) For bi-plane loops, Reference-Frontal is always selected by default at startup and after navigating from any loop to another. To select Reference-Lateral, click only in the Image Viewer area of the Lateral monitor.

A check appears in the active-plane indicator on the Reference monitor of the selected plane. For mono-plane studies (studies performed on a bi-plane system but only on one plane), an "X" appears in the active-plane indicator on the Reference monitor of the unused plane.

QCA is performed on the selected plane. Click anywhere in the viewing area of the desired Reference monitor to select that plane. On the side menu, click [Funct]-[QCA]. ISO center value is registered to calibration coefficient as an initial setting value. Manually perform the calibration if necessary, and then QCA analysis.

For LVA, click [Funct]-[LVA] on the side menu. Display the End Diastolic (ED) images and click anywhere in the viewing area of the desired Reference monitor on which you wish to perform calibration. ISO center value is registered to calibration coefficient as an initial setting value. Manually perform the calibration if necessary, and then LV analysis. The same calibration tool must be used for both planes. Each plane can be calibrated separately. Both calibration factors are displayed on the Reference monitors.

**Calibration Factor:**  
Fr 0.072 mm/pixel  
Lat 0.061 mm/pixel

Draw the RAO-ED contour on the monitor of one plane, and the LAO-ED contour on the monitor of the other plane. Then display the End Systolic (ES) image, and draw the RAO-ES contour on


the monitor of one plane, and the LAO-ES contour on the monitor of the other plane. Finally, click [Compute].

Image annotations can be made on one plane at a time. On the side menu, click [Funct]-[ABC Edit]. Click anywhere in the viewing area of the desired Reference monitor to select that plane. Make annotations on the selected plane. If you use the [Funct]-[Save] option to save your annotated image, only the image of the active plane is saved.

For Bi-plane loop, images are added to on both Frontal and Lateral plane at the same time when adding reference images. Reference images on each plane is displayed on each monitor. When the single-plane loop with a reference image is selected, other plane will stop. See ["4.6.2 Adding a Reference Image" P.4-60](#).

- **Switching to Independent Mode**

If it is desired to switch into Ref Independent mode, click [RefSync] button on the Reference monitor side menu (button no longer appears pushed in).



In Ref Independent mode, enable to perform actions independently on either the Frontal or Lateral Reference monitors. Image playback, playback speed, and zoom are performed only on selected plane monitor.

- **IVR NEO/SMART Touch**


The IVR NEO/SMART Touch [Monitor Select] button changes focus between the Reference and Acquisition monitors. For IVR NEO, when the button LED is on, the Acquisition monitors are selected. When it is off, the Reference monitors are selected. For SMART Touch, focused monitor is displayed ([ACQ] or [rEF]). Thereafter, commands are applied to the selected (focused) monitors.

If you select Ref Independent mode in the GUI, work with the IVR NEO on the selected plane.

### ■ Acquisition Monitor

- **Acq Synchronized Mode**

Click Acq Synchronized mode (AcqSync side-menu button pushed in) for Acq Synchronized mode. When this mode is selected, enable to operate both Frontal and Lateral Acquisition PC at the same time by pressing [Ctrl] key.



Press [Ctrl]+[->] or [Ctrl]+[-<]. The loop on the Frontal Acquisition monitor jumps to the next ([Ctrl]+[->]) or previous ([Ctrl]+[-<]) frame and pauses. The Lateral Acquisition monitor does the same (regardless of what frame it was showing before), pausing on the same frame number (but Lateral plane) as on the Frontal monitor.

Example:

1. Frontal is paused on frame 125 in loop 1, Lateral is paused on frame 75 in its loop 1.
2. Click [In sync] on the Reference monitor ide menu.
3. Press [Ctrl] and [->] on keyboard.
4. Frontal displays loop 1, frame 126, and Lateral displays its loop 1, frame 126.













And you can control playback speed by pressing [Ctrl] and mouse-clicking the speed control on the Reference monitor side menu. For other available features, see the next section "Acquisition Control with Keyboard and Mouse".



When in Acq Synchronized mode, [Ctrl] and [Ctrl]+[Alt] do the same thing. Press neither the [Ctrl] nor [Ctrl]+[Alt] keys to interact with the Reference monitor GUI.

- **Keyboard and Mouse**

Since there is no keyboard and mouse on both Acquisition computers, Acquisition GUI interaction is performed with the keyboard and mouse in one of two modes, "Acq Synchronized" or "Independent". The following table summarizes how the keyboard and mouse can be used to interact with the Acquisition monitor GUI when Frontal and Lateral are Synchronized. (If Independent, use [Ctrl] key for Frontal, and [Ctrl]+[Alt] keys for Lateral)


| No. | Action  | [Ctrl]+Keyboard or GUI   |
|-----|---|--|
| 1   | Play/show next loop/frame.  |   |
| 2   | Play/show previous loop/frame.  |   |
| 3   | If playing, pause.  |  or  |
| 4   | If paused, show next/previous frame.  |  /   |
| 5   | If paused, start playing. ( <i>Hold key for more than 2 seconds to begin forward / backward play.</i> ) |  /   |
| 6   | Playback speed: increase (right-click) / decrease (left-click).   |   |
| 7   | Brightness: increase (right-click) / decrease (left-click).   |   |
| 8   | Contrast: increase (right-click) / decrease (left-click).   |   |
| 9   | Zoom: increase (right-click) / decrease (left-click).   |   |



For No. 7 and 8 in the table above, disable to change the value of brightness and contrast on Frontal and Lateral at the same time by using keyboard though [AcqSync] button on the side menu is selected.

- **Switching to Acq Independent Mode**

In Acq Independent mode, use keyboard or mouse to perform actions independently on either the Frontal or Lateral Acquisition computers. Select the Frontal (A) Acquisition GUI by pressing the [Ctrl] key, and select the Lateral (B) Acquisition GUI by pressing the [Ctrl]+[Alt] keys together.



If it is desired to switch into Acq Independent mode, click [AcqSync] button on the Reference monitor side menu (button no longer appears pushed in).

- **IVR NEO/SMART Touch**

The IVR NEO/SMART Touch Monitor Select button changes focus between the Reference and Acquisition monitors. For IVR NEO, when the button LED is on, the Acquisition monitors are selected.

IVR NEO interaction is always synchronized regardless of Acq Synchronized/Independent mode settings in the GUI. If you select Acq Independent mode in the GUI and then work with the IVR NEO, it is switch to Acq Synchronized mode.




The Select and Split functions ([16-up Selector] and [Single/split Display] buttons respectively) of IVR NEO/SMART Touch are not available for the Acquisition-Lateral monitor.

### 10.2.4 Filtering Bi-plane Images by C-arm Position

To filter Image Selector icons for images acquired with a C-arm position that matches the current C-arm position (within a configurable tolerance), follow this procedure. Tolerance is configurable.

 "Software Tab" P.17-27.


- 1 Position the Frontal or Lateral C-arm, or both at the desired position.
- 2 Press the IVR NEO/SMART Touch [Filter by C-arm] button. The LED of IVR NEO turns on to indicate that Filter mode is selected.
- 3 Alternatively, click [Image Filtering] button on the bottom GUI bar and then in the Image Filtering list, choose [Match Fr+Lat C-Arm] mode.



The three filtering modes work as follows:


| Option                     | Description   |
|----------------------------|---|
| <b>Match Frontal C-Arm</b> | Reference-Frontal Image selector is filtered to show icons for loops that match Frontal C-arm position. The Lateral Image Selector is unaffected. |
| <b>Match Lateral C-Arm</b> | Reference-Lateral Image selector is filtered to show icons for loops that match Lateral C-arm position. The Frontal Image Selector is unaffected. |
| <b>Match FR+Lat C-Arm</b>  | Both Image Selectors are filtered to match the respective C-arm positions.  |

The Reference-Frontal and/or Reference-Lateral Image Selectors are filtered so only the icons of images that were acquired at the current C-arm position for the selected plane(s) are shown with the LAO/RAO and CRA/CAU angles at the top like this:




- 4 To revert to a non-filtered Image Selectors, press the IVR NEO/SMART Touch [Filter by C-Arm] button again (LED of IVR NEO turns off).

Or click the Image Filtering button on the bottom GUI bar and then, in the Image Filtering list, choose [All Images].




The Image Selectors are no-longer filtered; all icons are shown like this:



## 10.3 Bi-plane Image Preview

### Series Display on Studies Management Window

In the Series area of the Studies Management window, loops are identified as to the their plane, either **A** for Frontal or **B** for Lateral.




### GUI Interaction During Review

When no study is open for acquisition and reviewing a study on the Reference monitors, the [Ctrl] and [Ctrl-Alt] keys are not used. Reference monitor GUI interaction is performed with the keyboard and mouse. Loop playback is always synchronized between the Reference-Frontal and Reference-Lateral monitors. When a mono-plane loop is selected on the Reference-Frontal or Reference-Lateral monitor, the non-applicable plane loop is paused. Non-acquisition-related IVR NEO/SMART Touch functionality, including joystick playback control, is also available during review.



# 10.4 Shutdown

On the Reference monitor side menu, click [Exit]. All three computers are shutdown simultaneously. For additional details, see  ["3.2 User Authentication" P.3-4](#).

## 10.5 Administration

Several Options tabs have additional values related to the Bi-plane option.

Bi-plane options are configured on the two Acquisition computer [System] tabs as follows:

- On the Frontal [System] tab, in the Link box, Biplane and (A) Frontal are selected and the Computer Name of the Lateral Acquisition computer is entered in Plane B Name.

The screenshot shows the 'System Configuration' dialog box. On the left is a sidebar with various system settings. The main area is divided into 'System Information' and 'Link' sections. In the 'Link' section, 'Link Mode' is set to 'Biplane' (checked), and '(A) Frontal' is also selected (radio button). The 'Plane B Name' field contains the text 'DAR9500ACQ-L-01'. Other fields in 'System Information' include Computer Name, Institution Name, Station Name, Machine Serial Number, DICOM AE Title, Manufacturer Model Name, Implementation UID, and Information Version.

- On the Lateral [System] tab, Biplane and (B) Lateral are selected.

This screenshot shows the 'System Configuration' dialog box for the Lateral Acquisition computer. In the 'Link' section, 'Link Mode' is set to 'Biplane' (checked), and '(B) Lateral' is selected (radio button). The 'Plane B Name' field is empty. The 'System Information' section contains similar fields to the previous screenshot, with some values updated to reflect the Lateral configuration.

- On DUP setting window of the Reference monitor, select [Biplane] in [Program Description] to set Bi-plane program.

The screenshot shows the 'Program Description' window. It has a text field containing 'DA' and four checkboxes: 'Frontal', 'Lateral', 'Biplane', and 'Timer'. The 'Biplane' checkbox is circled in red, indicating it should be selected.

See also ["17.7.4 Menus and DUP Configuration" P.17-69.](#)

- On the Reference monitor [Software] tab, Bi-plane options are configured as follows:

**Configuration**

**Software Configuration**

**Close Study**

Set default to ☒ Close study in Acquisition ☒ Close study in Review

**Print Options**

☒ Print in Color (Non-DICOM) ☒ Keep the Print Manager open after printing

**Image Selector Filtering (from C-Arm position)**

Degree of tolerance

**View Direction**

Frontal

Lateral

**Automatic Plane Mode Change**

☒ Enable

**Remote Maintenance**


☒ Site-View Plus

URL

☒ MPC Error Report

MPC

OK Cancel Apply

| Item                               | Description  | Modify |
|------------------------------------|--|--------|
| <b>View Direction</b>              | Frontal: Set Frontal eye direction.<br>Lateral: Set Lateral eye direction.   | ×      |
| <b>Automatic Plane Mode Change</b> | <p><b>Automatically change the plane depending on the C-arm condition. Ex.)</b></p> <ul style="list-style-type: none"> <li>When changing the C-arm mode from Bi-plane to Single-plane, it will be changed to Frontal or Lateral.</li> <li>When changing the C-arm mode from Single-plane to Bi-plane, acquisition plane will be changed to "Bi".</li> </ul> <p><b>Do not check this item usually.</b></p> <div style="border: 1px solid black; padding: 10px; margin-top: 10px;"> <p style="text-align: center;"> <b>NOTE</b></p> <p>When changing the C-arm mode to Bi-plane, acquisition plane will surely be changed to "Bi". Do not valid this setting if there are many cases to acquire images on Frontal or Lateral after changing to "Bi".</p> </div> | ○      |

## 10 Difference with Bi-plane Option

- On the Reference monitor [Hardware] and [Hardware (Shutter)] tab, Bi-plane options are configured as follows:

**Configuration**

**Hardware Configuration**

System  
DICOM Hosts  
DICOM  
Software  
**Hardware**  
Hardware (Shutter)  
Storage (Local)  
Storage (Network)  
Display  
Menus  
Notifications  
Physicians  
Study Information  
Fonts & Colors  
Database  
Devices  
External Software  
Fusion  
Logs

**Information**

|                         | Frontal | Lateral |
|-------------------------|---------|---------|
| Grid                    | IN      | IN      |
| FPD Size (inches)       | 12      | 12      |
| FOD (cm)                | 72      | 72      |
| ISO Center To Skin (cm) | 15      | 15      |

**Field Of View (inches)**

|          | Frontal | Lateral | FOV 5th | FOV 6th | FOV 7th | FOV 8th |
|----------|---------|---------|---------|---------|---------|---------|
| Large    | 12      | 12      | 4.5     | 4.5     |         |         |
| Medium   | 10      | 10      | 0       | 0       |         |         |
| Small    | 8       | 8       | 0       | 0       |         |         |
| Smallest | 6       | 6       | 0       | 0       |         |         |

**Dosemeter Configuration**

Dosemeter: VacuDAQ

**Table Configuration**

☒ Table Reposition

☒ Show ROI indicating the table position

Stroke Length (mm): 1350

**Large Monitor Configuration**

☒ Large Monitor Control

LMM IP: 192.168.100.37

**Curves Acquisition**

| Channel                                       | Label    |
|---|----------|
| <input checked="" type="checkbox"/> Channel 1 | ECG      |
| <input type="checkbox"/> Channel 2            | Pressure |
| <input type="checkbox"/> Channel 3            | PHYSIO   |
| <input type="checkbox"/> Channel 4            | FLOW     |

☒ Default Acquire Waves

**SMART Touch Configuration**

|   | Unit        |
|---|-------------|
| <input checked="" type="checkbox"/> Touch 1 | TOUCH-01-01 |
| <input checked="" type="checkbox"/> Touch 2 | TOUCH-01-02 |
| <input checked="" type="checkbox"/> Touch 3 | TOUCH-01-03 |

OK Cancel Apply

**Configuration**

**Hardware (Shutter) Configuration**

System  
DICOM Hosts  
DICOM  
Software  
Hardware  
**Hardware (Shutter)**  
Storage (Local)  
Storage (Network)  
Display  
Menus  
Notifications  
Physicians  
Study Information  
Fonts & Colors  
Database  
Devices  
External Software  
Fusion  
Logs


**Shutter Positions**



|          | Frontal   |       |           |       |           |       | Lateral   |       |           |       |           |       |
|----------|-----------|-------|-----------|-------|-----------|-------|-----------|-------|-----------|-------|-----------|-------|
|          | C-Shutter |       | H-Shutter |       | V-Shutter |       | C-Shutter |       | H-Shutter |       | V-Shutter |       |
|          | Open      | Close | Open      | Close | Open      | Close | Open      | Close | Open      | Close | Open      | Close |
| Large    | -15       | 62    | -15       | 100   | -15       | 100   | -15       | 62    | -15       | 100   | -15       | 100   |
| Medium   | -16       | 56    | -15       | 100   | -15       | 100   | -16       | 56    | -15       | 100   | -15       | 100   |
| Small    | -18       | 50    | -15       | 100   | -15       | 100   | -18       | 50    | -15       | 100   | -15       | 100   |
| Smallest | -24       | 35    | -15       | 100   | -15       | 100   | -24       | 35    | -15       | 100   | -15       | 100   |
| FOV 5th  | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   |
| FOV 6th  | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   |
| FOV 7th  | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   |
| FOV 8th  | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   | 0         | 100   |

OK Cancel Apply

| Item                 | Description  | Modify |
|----------------------|--|--------|
| <b>Grid</b>          | Defines whether a grid is used in the Frontal and Lateral X-ray imaging system. Select <b>IN</b> if a grid is used. This option is used to complete the DICOM information for acquired images. This option must match the installation.  | ×      |
| <b>FOD (cm)</b>      | Specifies the focal point to object (ISO Center) distance for the Frontal and Lateral planes. Also known as SOD (Source to Object Distance). This is a fixed value (typically 72cm) according to C-arm model. This value (cm) must match the installation. (0 to 100.)   | ×      |
| <b>Field of View</b> | Defines eight field of view values (inches) for each plane. This option is used to complete the DICOM information for acquired images. This option must match the installation. <b>Large</b> can be no larger than <b>FPD Size</b> . <b>Medium</b> , <b>Small</b> , and <b>Smallest</b> must be progressively smaller. | ×      |

- Settings of [Display] tab for Bi-plane option are shown as follows.



| Item               | Description   | Modify                 |
|--------------------|---|------------------------|
| Virtual Collimator | <p>Virtual collimators are displayed on the Acquisition monitor showing the position of the hardware collimators.</p> <p><b>Delay (sec)</b> defines the time in seconds (0 to 100) that virtual collimators are displayed after the hardware collimators have stopped moving.</p> <p>Width defines the thickness of the displayed line (range 2 to 30).</p> <div> <b>NOTE</b></div> <p>Virtual collimation function indicates the collimator leaf position before fluoroscopy as a suggestion. Therefore, the position may misalign by the field of view size, C-arm position, and etc.</p> <div> <b>NOTE</b></div> <p>[For Bi-plane] When operating Virtual collimator from the local console in examination room, Virtual collimator will be displayed on the acquisition monitor which performed the last fluoroscopy.</p> | <div><div></div></div> |

See also  "17.7 Options Configuration" P.17-14.

This page is intentionally left blank.

# *Chapter 11*

---

## SCORE StentView (Option)

### ■ Description

|                            |       |
|----------------------------|-------|
| 11.1 Overview .....        | 11-2  |
| 11.2 Specification .....   | 11-3  |
| 11.3 Handling Method ..... | 11-4  |
| 11.4 Troubleshooting ..... | 11-18 |

# 11.1 Overview

SCORE StentView Option (hereafter called “StentView”) is an attachment function to DAR-9500f SHIMADZU digital angiography system. Extract the marker position from real-time image and highlight the device by digital image processing. When using two Stent, this function is good for confirming the second Stent position.

## 11.1.1 Features

---

StentView has the following features.

- Apply the processing to real-time acquisition image and display the device highlighted image in live image.
- Enable to install on DAR-9500f standard system. Work station is not necessary.



# 11.2 Specification

## 11.2.1 Specification of SCORE StentView

|                  |   |
|------------------|---|
| Acquisition Rate | 15fps, 30fps  |
| Image Size       | <Full Screen View><br>(Whole)<br>1024x1024 matrix   |
| Image Processing | <ul style="list-style-type: none"> <li>• Automatic gray level adjustment</li> <li>• Edge enhancement processing</li> <li>• ROI setting</li> </ul> |

## 11.2.2 Features

### 11.2.2.1 Image Acquisition

| No. | Items     | Features  |
|-----|-----------|---|
| 1   | StentView | Extract a marker position from a real-time image and enhance and display the device by digital image processing. Add some frames and display enhanced stent image.<br>This is useful when place 2 stent and confirm the position of second one. |

### 11.2.2.2 Image Processing

| No. | Items                      | Features  |
|-----|----------------------------|---|
| 1   | Auto Window Level          | Images are displayed with automatically-stabilized contrast.  |
| 2   | Real-time Edge Enhancement | Sharp images are acquired by emphasizing the subject edge with spatial frequency emphasis.  |
| 3   | ROI Setting                | If images have more than three points of markers or include high-contrast artificial structures, enable to detect markers correctly and display the result by using ROI setting to remove unnecessary area. |

## 11.3 Handling Method

### 11.3.1 Introduction

---


By pressing StentView button to enable StentView processing for the next radiography. StentView button is enabled only when StentView is selected on Digital User Program (hereafter called “DUP”).

### 11.3.2 Displaying StentView

---

#### ■ Full Screen View

During StentView processing, a zoomed StentView image on the ACQ monitor and X-ray exposed whole image on the REF monitor will be displayed.



ACQ Monitor




REF Monitor

### 11.3.3 StentView Processing

---

Follow the instructions below to apply StentView processing.


- 1 Confirm the current radiography digital program setting displayed at the bottom of image viewer window. Be sure that program is set to StentView.



- 2 Position the marker. Change the C-arm angle, bed position and FOV size to position the marker to be fit in the image.

- 3 Click [Program: StentView] on the side menu or press [StentView] button on IVR NEO/SMART Touch for StentView radiography.

\*If StentView button is assigned on Keyboard shortcut, enable StentView radiography from the shortcut key.



Side Menu



IVR NEO



SMART Touch



When selecting unassigned StentView option, RAD StentView button cannot be used.



- 4 Press the radiography switch (foot switch or hand switch). StentView radiography is started and StentView image is displayed.



## WARNING

If unable to detect markers normally during StentView radiography, display the last image detected markers and "LIH" will display at the bottom-right of ACQ monitor. StentView result is wrong if markers are not detected normally. Stop using StentView option immediately.

Refer to "[11.4.1 Phenomenon and Action](#)" for handling problems.



## WARNING

Two marker points are required for recognition for StentView. If a point and more than three marker points are existed within ROI, it gives not correct results.



## WARNING

If contrasty artificial structure such as a pacemaker and a prosthetic valve is placed within the same ROI as a marker, structure may fail to recognize as a marker and may not obtain correct result.



## NOTE

As displaying StentView by adding a number of frames, it cannot obtain correct results if the position of marker and device, and X-ray exposure angle has changed. And if contrast media is injected, the false results may be displayed because the marker is hidden with the contrast media.

- 5 After the image acquisition, release hand switch or foot switch.  
When releasing a switch, automatically a marker detected frame starts playing.

- 6 If you want to continue StentView radiography, go back to step 2.

When completing StentView radiography, it returns to the normal radiography program automatically. If you want to continue StentView radiography, repeat procedures from step 2.
















Measurement tools such as QCA (Quantitative Coronary Artery Analysis) and LV (Left Ventricle Analysis) cannot use for StentView radiography acquisition image.

### 11.3.3.1 Setting ROI

If marker does not detect normally because of body or artificial structures similar to marker shape, set one or two ROI to reduce the possibility of false detection. Follow the instructions below to set ROI.

- 1 Perform fluoroscopy at least two seconds and check the change of marker position by heartbeat.
- 2 Click ROI setting button of StentView on the side menu or IVR NEO/SMART Touch and ROI is displayed on ACQ monitor which is currently selected as an acquisition plane. When ROI is displayed on LIH (Last Image Hold), automatically replays the last fluoroscopy image.

|   |   |   |   |  |
|---|---|---|---|--|
| GUI                                     |    |   |   |  |
|   | <b>Button</b>   |   | <b>Description</b>  |  |
|   |    |   | Displays when ROI is not setting.<br>Click to display one ROI.                                  |  |
|   |    |   | Displays when one ROI is displayed.<br>Click to display two ROI.                                |  |
|   |   |   | Displays when two ROI are displayed.<br>Click to disable ROI and enables normal StentView mode. |  |
| IVR NEO/<br>SMART<br>Touch/<br>Keyboard |  |  |              | Press a button to change the number of ROI as 1-2-0.                           |
|   |  |  |              | Press a button to display one ROI. If already displayed, ROI will be disabled. |
|   |  |  |              | Press a button to display two ROI. If already displayed, ROI will be disabled. |

- 3 Set ROI on the ACQ monitor. Adjust the number, position and size of ROI if necessary.




ROI must include all range of movement of marker by heartbeat to get correct result. Set the number, position and size of ROI as checking fluoroscopy image replaying on the Acquisition monitor.

- **Setting Numbers**


- **[1-ROI]**

- Set ROI including the range of movement of two markers.



- **[2-ROI]**

- Set two ROI for Proximal and Distal markers. It is no problem if those two ROI overlap, but one marker must be included.



- **Setting Position**

- [Mouse]**

- Move the mouse pointer to the ACQ monitor and click, drag and drop within the ROI to modify the position.

- [IVR NEO/SMART Touch]**

- ROI does activate with move mode. When tilting the joystick of IVR NEO up/down or right/left, ROI moves to the tilting direction.

- When pressing the joystick button, ROI mode is changed to size modification mode. Each mode will change as follows.

- [1-ROI]**

- Move mode-->Size Modification mode-->Move mode

- [2-ROI]**

- Proximal Move mode-->Proximal Size Modification mode-->Distal Move mode-->Distal Size Modification mode-->Proximal Move mode.

Enable to distinguish the move mode and size modification mode by shape of ROI.

To change the Lateral position, press [Lateral Operation] button and then set the position.

- **Setting Size**

- [Mouse]**

- Move the mouse pointer to the ACQ monitor and click, drag and drop within the ROI to modify the size.

- [IVR NEO/SMART Touch]**

- When ROI is in the move mode, ROI mode is changed to size modification mode by pressing the joystick button.

- During size modification mode, tilt the joystick to up/down and right/left to fix the point at the top-left of ROI, and the point at the bottom-right of ROI does activate to change the size.

To change the Lateral size, press [Lateral Operation] button and then set the size.

- 4 Press the radiography switch (foot or hand switch). StentView radiography starts and StentView image will display.




11.3.4 Change View

On ACQ and REF monitor, enable to change the display from StentView to original image by canceling a special image processing StentView.

Click [Change] button on the side menu to change the image. It is start playing from the 1st frame when playing an original StentView image.

\*If change StentView is assigned to keyboard shortcut and IVR NEO/SMART Touch button, enable to change StentView with a shortcut key.



11.3.5 Setting StentView




Configure StentView on REF PC.

## Setting Shortcut Keys


Set IVR NEO and Shortcut keys with [Shortcut keys]/[SMART Touch] on [Physicians] tab. Only Shimadzu or specified service representatives can change the settings.

Refer to ["17.7 Options Configuration"](#) for setting procedure.

Selecting Shortcut keys:



Selecting SMART Touch:




Description of StentView/StentShot of [Shortcut keys] and [SMART Touch] tab is shown as follows.

| No.                       | Item                               | Description   |
|---------------------------|------------------------------------|---|
| Shortcut keys/SMART Touch |                                    |   |
| StentView/StentShot       |                                    |   |
| 1                         | Activate StentView/StentShot       | Change Enable/Disable of StentView/StentShot.   |
| 2                         | Activate ROI-StentView/StentShot   | Change disable of 1-ROI StentView/2-ROI StentView<br>1-ROI StentShot/2-ROI StentShot. |
| 3                         | Activate 1-ROI StentView/StentShot | Change Enable/Disable of 1-ROI StentView/StentShot.                                   |
| 4                         | Activate 2-ROI StentView/StentShot | Change Enable/Disable of 2-ROI StentView/StentShot.                                   |
| 5                         | Change View                        | Change StentView/StentShot image view/Original image view.                            |

## Setting StentView

Set StentView/StentShot with [PCI App.] on [Physicians] tab.




Description of [PCI App.] tab is shown as follows.

| No.                 | Item                            | Description   |
|---------------------|---------------------------------|---|
| StentView/StentShot |                                 |   |
| 1                   | Coordination with Geometry      | Check this item to disable ROI when changing geometry.  |
| 2                   | 1-ROI Initial Size (%)          | During study, enable to register the initial size of first 1-ROI.<br>Enable to set value at rate of 100% of the whole image as 100%.  |
| 3                   | 2-Proximal ROI Initial Size (%) | During study, enable to register the Proximal initial size of first 2-ROI.  |
| 4                   | 2-Distal ROI Initial Size (%)   | During study, enable to register the Distal initial size of first 2-ROI.  |
| 5                   | ROI Minimum Size (%)            | Enable to set minimum size of ROI.  |
| 6                   | Hold the last ROI status        | Check this item to save the last ROI position and size, and enable to use the same ROI for next StentView.<br>However, if mechanical equipment such as tabletop and C-arm operate, ROI position and size will be initialized even though this item is checked.<br>Uncheck this item to display ROI always in the center of the image. |
| 7                   | Auto-Replay last Fluoro         | Check this item to replay the last fluoro for 2 seconds automatically, when ROI is displayed on LIH.  |
| 8                   | Scale of StentView area         | Enable to change the magnification percentage of StentView area.<br>Normally, do not change the value which is 1.3 times of default percentage.   |

| No. | Item   | Description   |
|-----|--|---|
| 9   | Scale of full screen StentView/<br>StentShot | Enable to change the magnification percentage of full screen StentView/StentShot. Normally, do not change the value which is 1.5 times of default percentage. |
| 10  | FullScreen Mode                              | Check this item to display StentView area in full screen.   |

## Geometry

Enable to set margins that automatically disable ROI with [Geometry] on [Physicians] tab.



Description of [Geometry] tab is shown as follows.

| No.                        | Item   | Description  |
|----------------------------|--|--|
| Coordination with Geometry |  |  |
| 1                          | Guide Display Timeout (sec.)                     | Not used.  |
| 2                          | Oblique Angle Margin (degree)                    | If Coordination with Geometry is checked and the C-arm is changed more than input angle in an oblique direction, SIMAP and Guide Line will hide. |
| 3                          | Sagittal Angle Margin (degree)                   | If Coordination with Geometry is checked and the C-arm is changed more than input angle in a sagittal direction, SIMAP and Guide Line will hide. |
| 4                          | SID Margin (mm)                                  | If Coordination with Geometry is checked and SID is changed more than input value, SIMAP and Guide Line will hide.                               |
| 5                          | Table Longitudinal Position Margin (mm)          | If Coordination with Geometry is checked and the table moves more than input value in a longitudinal direction, SIMAP and Guide Line will hide.  |
| 6                          | Table Transversal Position Margin                | If Coordination with Geometry is checked and the table moves more than input value in a transversal direction, SIMAP and Guide Line will hide.   |
| 7                          | Table Height Margin                              | If Coordination with Geometry is checked and the table height is changed more than input value, SIMAP and Guide Line will hide.                  |
| 8                          | Ceiling Travel Longitudinal Position Margin (mm) | If Coordination with Geometry is checked and the C-arm moves more than input value in a longitudinal direction, SIMAP and Guide Line will hide.  |

| No. | Item  | Description  |
|-----|---|--|
| 9   | Ceiling Travel Transversal Position Margin (mm) | If Coordination with Geometry is checked and the C-arm moves more than input value in a transversal direction, SIMAP and Guide Line will hide. |

# 11.4 Troubleshooting

## 11.4.1 Phenomenon and Action

When StentView does not work normally, perform following actions.

| Phenomenon  |
|---|
| Nothing appears on StentView area.  |
| Assumed Cause   |
| Disable to set marker position.   |
| Action  |
| <p>If disabling to set marker, there are some possible cases. Perform following actions correspond to the situation.</p> <p>[Case 1] Marker is placed outside the window because of heartbeat.<br/>Change FOV size and bed position to fit the marker within the display area.</p> <p>[Case 2] Fast heartbeat.<br/>The faster the marker moves by heartbeat, the more the marker detection rate goes down.<br/>Perform StentView radiography with higher exposure rate.</p> <p>[Case 3] Large FOV size.<br/>The larger the FOV size gets, the smaller the marker size gets and detection rate goes down.<br/>Set FOV size smaller or set ROI within the range of movement of the marker.</p> <p>[Case 4] Large C-arm angle.<br/>The deeper the C-arm angle gets, the more the marker contrast and detection rate goes down.<br/>Reach performance limits. Perform StentView radiography with shallower angle.</p> |

| Phenomenon  |
|---|
| <p>When showing StentView image, delay a showing of next/previous image.</p> <p>Slow reaction and operation of a mouse.</p> <p>Image does not transfer from ACQ to REF.</p> |
| Assumed Cause   |
| PC is highly-load excessively.  |
| Action  |
| <p>Do not perform StentView radiography for a while (a minute or two).</p> <p>If it does not improve for a long time, restart the system.</p>                               |



|  |
|--|
| Phenomenon   |
| StentView does not activate, even when pressing a hand switch.   |
| Assumed Cause  |
| Software error, hardware malfunction.  |
| Action   |
| Restart the system. If it is not improved after restarting the system, perform normal fluoroscopy and radiography, and contact our service agency. |

|  |
|--|
| Phenomenon   |
| Abnormal image of StentView radiography is shown.  |
| Assumed Cause  |
| StentView parameter is not setting correctly.  |
| Action   |
| Select other DUP once and then select StentView DUP again.<br>If abnormal image is shown for normal fluoroscopy, restart the system. |

|   |
|---|
| Phenomenon  |
| Result does not display correctly on the StentView area though setting ROI. |
| Assumed Cause   |
| Marker is placed outside the ROI.   |
| Action  |
| Set ROI including the range of movement of marker by heartbeat.             |

This page is intentionally left blank.

# Chapter 12

---

## SCORE StentShot (Option)

### ■ Description

|                            |       |
|----------------------------|-------|
| 12.1 Overview .....        | 12-2  |
| 12.2 Specification .....   | 12-3  |
| 12.3 Handling Method ..... | 12-4  |
| 12.4 Troubleshooting ..... | 12-18 |

## 12.1 Overview

SCORE StentShot Option (hereafter called “StentShot”) is an attachment function to DAR-9500f SHIMADZU digital angiography system. Extract the marker position from real-time image and highlight the device by digital image processing. When using two Stent, this function is good for confirming the shape of placed Stent.

### 12.1.1 Features

---

Stent has the following features.

- Apply the processing to real-time acquisition image and display the device highlighted image in live image.
- Enable to install on DAR-9500f standard system. Work station is not necessary.

# 12.2 Specification

## 12.2.1 Specification of SCORE StentShot

|                  |   |
|------------------|---|
| Acquisition Rate | 15fps, 30fps  |
| Image Size       | <Full Screen View><br>(Whole)<br>1024x1024 matrix   |
| Image Processing | <ul style="list-style-type: none"> <li>• Automatic gray level adjustment</li> <li>• Edge enhancement processing</li> <li>• ROI setting</li> </ul> |

## 12.2.2 Features

### 12.2.2.1 Image Acquisition

| No. | Items     | Features   |
|-----|-----------|--|
| 1   | StentShot | Extract a marker position from a real-time image, and enhance and display the device by digital image processing. Add some frames and display enhanced stent images.<br>This is useful for confirming the shape of placed stent. |

### 12.2.2.2 Image Processing

| No. | Items                      | Features  |
|-----|----------------------------|---|
| 1   | Auto Window Level          | Images are displayed with automatically-stabilized contrast.  |
| 2   | Real-time Edge Enhancement | Sharp images are acquired by emphasizing the subject edge with spatial frequency emphasis.  |
| 3   | ROI Setting                | If images have more than three points of markers or include high-contrast artificial structures, enable to detect markers correctly and display the result by using ROI setting to remove unnecessary area. |

## 12.3 Handling Method

### 12.3.1 Introduction

---


By pressing StentShot button to enable StentShot processing for the next radiography. StentShot button is enabled only when StentShot is selected on Digital User Program (hereafter called “DUP”).

### 12.3.2 Displaying StentShot

---

#### ■ Full Screen View

During StentShot processing, display zoom in StentShot image on ACQ monitor and X-ray exposed whole image on REF monitor.



ACQ Monitor



REF Monitor

### 12.3.3 StentShot Processing

---


Follow the instructions below to apply StentShot processing.



As displaying StentShot by adding a number of frames, it cannot obtain correct results if the position of marker and device has changed between frames.


Do not operate catheter during StentShot.

- 1 Confirm the current radiography digital program setting displayed at the bottom of image viewer window. Be sure that program is set to StentShot.




- 2 Position the marker. Change the C-arm angle, bed position and FOV size to position the marker to be fit in the image.
- 3 Click [Program: StentShot] on the side menu or press [StentShot] button on IVR NEO/SMART Touch for StentShot radiography.

\*If StentShot button is assigned on Keyboard shortcut, enable StentShot radiography from the shortcut key.



When selecting unassigned StentShot option, RAD StentShot button cannot be used.



- 4 Press the radiography switch (foot switch or hand switch). StentShot radiography is started and StentShot image is displayed.



### WARNING

If unable to detect markers normally during StentShot radiography, display the last image detected markers and "Cannot detect target devices" will display at the top of ACQ monitor. StentShot result is wrong if markers are not detected normally. Stop using StentShot option immediately.

Refer to ["12.4.1 Phenomenon and Action"](#) for handling problems.



### WARNING

Two marker points are required for recognition for StentShot. If a point and more than three marker points are existed within ROI, it gives not correct results.



### WARNING

If contrasty artificial structure such as a pacemaker and a prosthetic valve is placed within the same ROI as a marker, structure may fail to recognize as a marker and may not obtain correct result.



### NOTE

As displaying StentShot by adding a number of frames, it cannot obtain correct results if the position of marker and device, and X-ray exposure angle has changed. And if contrast media is injected, the false results may be displayed because the marker is hidden with the contrast media.



## 5 After the image acquisition, release hand switch or foot switch.

When releasing a switch, the last frame will be displayed. For an original image, transfer the image to the image server as usual radiography image. And to send an image-processed image, only the last frame will be sent.



After StentShot radiography, do not start image acquisition until display the last frame on ACQ monitor

## 6 If you want to continue StentShot radiography, go back to step 2.

When completing StentShot radiography, it returns to the normal radiography program automatically. If you want to continue StentShot radiography, repeat procedures from step 2.



Measurement tools such as QCA (Quantitative Coronary Artery Analysis) and LV (Left Ventricle Analysis) cannot use for StentShot radiography acquisition image.

### 12.3.3.1 Setting ROI














If marker does not detect normally because of body or artificial structures similar to marker shape, set one or two ROI to reduce the possibility of false detection. Follow the instructions below to set ROI.

#### 1 Perform fluoroscopy at least two seconds and check the change of marker position by heartbeat.

#### 2 Click ROI setting button of StentShot on the side menu or IVR NEO/SMART Touch and ROI is displayed on ACQ monitor which is currently selected as an acquisition plane.

When ROI is displayed on LIH (Last Image Hold), automatically replays the last fluoroscopy image.

When ROI is displayed on StentShot image, automatically replays the original Image.

|   |   |   |   |  |
|---|---|---|---|--|
| GUI                                     |    |   |   |  |
|   | <b>Button</b>   |   | <b>Description</b>  |  |
|   |    |   | Displays when ROI is not setting.<br>Click to display one ROI.                                  |  |
|   |    |   | Displays when one ROI is displayed.<br>Click to display two ROI.                                |  |
|   |  |   | Displays when two ROI are displayed.<br>Click to disable ROI and enables normal StentShot mode. |  |
| IVR NEO/<br>SMART<br>Touch/<br>Keyboard |  |  |              | Press a button to change the number of ROI as 1-2-0.                           |
|   |  |  |              | Press a button to display one ROI. If already displayed, ROI will be disabled. |
|   |  |  |              | Press a button to display two ROI. If already displayed, ROI will be disabled. |

- 3 Set ROI on the ACQ monitor. Adjust the number, position and size of ROI if necessary.




ROI must include all range of movement of marker by heartbeat to get correct result. Set the number, position and size of ROI as checking fluoroscopy image replaying on the Acquisition monitor.

- **Setting Numbers**


- **[1-ROI]**

- Set ROI including the range of movement of two markers.



- **[2-ROI]**

- Set two ROI for Proximal and Distal markers. It is no problem if those two ROI overlap, but one marker must be included.



- **Setting Position**

- [Mouse]**

- Move the mouse pointer to the ACQ monitor and click, drag and drop within the ROI to modify the position.

- [IVR NEO/SMART Touch]**

- ROI does activate with move mode. When tilting the joystick of IVR NEO up/down or right/left, ROI moves to the tilting direction.

- When pressing the joystick button, ROI mode is changed to size modification mode. Each mode will change as follows.

- [1-ROI]**

- Move mode-->Size Modification mode-->Move mode

- [2-ROI]**

- Proximal Move mode-->Proximal Size Modification mode-->Distal Move mode-->Distal Size Modification mode-->Proximal Move mode.

Enable to distinguish the move mode and size modification mode by shape of ROI.

To change the Lateral position, press [Lateral Operation] button and then set the position.

- **Setting Size**

- [Mouse]**

- Move the mouse pointer to the ACQ monitor and click, drag and drop within the ROI to modify the size.

- [IVR NEO/SMART Touch]**

- When ROI is in the move mode, ROI mode is changed to size modification mode by pressing the joystick button.

- During size modification mode, tilt the joystick to up/down and right/left to fix the point at the top-left of ROI, and the point at the bottom-right of ROI does activate to change the size.

To change the Lateral size, press [Lateral Operation] button and then set the size.


- 4** Press the radiography switch (foot or hand switch). StentShot radiography starts and StentShot image will display.

12.3.4 Change View

On ACQ and REF monitor, enable to change the display from StentShot to original image by canceling a special image processing StentShot.

Click [Change] button on the side menu to change the image. It is start playing from the 1st frame when playing an original StentShot image.

\*If change StentShot is assigned to keyboard shortcut and IVR NEO/SMART Touch button, enable to change StentShot with a shortcut key.



12.3.5 Setting StentShot




Configure StentShot on REF PC.

## Setting Shortcut Keys


Set IVR NEO and Shortcut keys with [Shortcut keys]/[SMART Touch] on [Physicians] tab. Only Shimadzu or specified service representatives can change the settings.

Refer to ["17.7 Options Configuration"](#) for setting procedure.

Selecting Shortcut keys:



Selecting SMART Touch:




Description of StentView/StentShot of [Shortcut keys] and [SMART Touch] tab is shown as follows.

| No.                       | Item                               | Description   |
|---------------------------|------------------------------------|---|
| Shortcut keys/SMART Touch |                                    |   |
| StentView/StentShot       |                                    |   |
| 1                         | Activate StentView/StentShot       | Change Enable/Disable of StentView/StentShot.   |
| 2                         | Activate ROI-StentView/StentShot   | Change disable of 1-ROI StentView/2-ROI StentView<br>1-ROI StentShot/2-ROI StentShot. |
| 3                         | Activate 1-ROI StentView/StentShot | Change Enable/Disable of 1-ROI StentView/StentShot.                                   |
| 4                         | Activate 2-ROI StentView/StentShot | Change Enable/Disable of 2-ROI StentView/StentShot.                                   |
| 5                         | Change View                        | Change StentView/StentShot image view/Original image view.                            |

## Setting StentShot

Set StentView/StentShot with [StentView/StentShot] on [Physicians] tab.



Description of [StentView/StentShot] tab is shown as follows.


| No.                 | Item                            | Description   |
|---------------------|---------------------------------|---|
| StentView/StentShot |                                 |   |
| 1                   | Coordination with Geometry      | Check this item to disable ROI when changing geometry.  |
| 2                   | 1-ROI Initial Size (%)          | During study, enable to register the initial size of first 1-ROI.<br>Enable to set value at rate of 100% of the whole image as 100%.  |
| 3                   | 2-Proximal ROI Initial Size (%) | During study, enable to register the Proximal initial size of first 2-ROI.  |
| 4                   | 2-Distal ROI Initial Size (%)   | During study, enable to register the Distal initial size of first 2-ROI.  |
| 5                   | ROI Minimum Size (%)            | Enable to set minimum size of ROI.  |
| 6                   | Hold the last ROI status        | Check this item to save the last ROI position and size, and enable to use the same ROI for next StentShot.<br>However, if mechanical equipment such as tabletop and C-arm operate, ROI position and size will be initialized even though this item is checked.<br>Uncheck this item to display ROI always in the center of the image. |
| 7                   | Auto-Replay last Fluoro         | Check this item to replay the last fluoro for 2 seconds automatically, when ROI is displayed on LIH.  |
| 8                   | Scale of StentView area         | Enable to change the magnification percentage of StentView area.<br>Normally, do not change the value which is 1.3 times of default percentage.   |



| No. | Item   | Description   |
|-----|--|---|
| 9   | Scale of full screen StentView/<br>StentShot | Enable to change the magnification percentage of full screen.<br>Normally, do not change the value which is 1.5 times of default<br>percentage. |
| 10  | FullScreen Mode                              | Check this item to display StentShot area in full screen.   |

## Geometry

Enable to set margins that automatically disable ROI with [Geometry] on [Physicians] tab.



Description of [Geometry] tab is shown as follows.

| No.                        | Item   | Description  |
|----------------------------|--|--|
| Coordination with Geometry |  |  |
| 1                          | Guide Display Timeout (sec.)                     | Not used.  |
| 2                          | Oblique Angle Margin (degree)                    | If Coordination with Geometry is checked and the C-arm is changed more than input angle in an oblique direction, SIMAP and Guide Line will hide. |
| 3                          | Sagittal Angle Margin (degree)                   | If Coordination with Geometry is checked and the C-arm is changed more than input angle in a sagittal direction, SIMAP and Guide Line will hide. |
| 4                          | SID Margin (mm)                                  | If Coordination with Geometry is checked and SID is changed more than input value, SIMAP and Guide Line will hide.                               |
| 5                          | Table Longitudinal Position Margin (mm)          | If Coordination with Geometry is checked and the table moves more than input value in a longitudinal direction, SIMAP and Guide Line will hide.  |
| 6                          | Table Transversal Position Margin                | If Coordination with Geometry is checked and the table moves more than input value in a transversal direction, SIMAP and Guide Line will hide.   |
| 7                          | Table Height Margin                              | If Coordination with Geometry is checked and the table height is changed more than input value, SIMAP and Guide Line will hide.                  |
| 8                          | Ceiling Travel Longitudinal Position Margin (mm) | If Coordination with Geometry is checked and the C-arm moves more than input value in a longitudinal direction, SIMAP and Guide Line will hide.  |

| No. | Item  | Description  |
|-----|---|--|
| 9   | Ceiling Travel Transversal Position Margin (mm) | If Coordination with Geometry is checked and the C-arm moves more than input value in a transversal direction, SIMAP and Guide Line will hide. |

## 12.4 Troubleshooting

### 12.4.1 Phenomenon and Action

When StentShot does not work normally, perform following actions.

| Phenomenon  |
|---|
| Normal image displays after the 9th frame.<br>Does not update StentShot image and "Cannot detect target devices" is displayed at the top of ACQ monitor.  |
| Assumed Cause   |
| Disable to set marker position.   |
| Action  |
| If disabling to set marker, there are some possible cases. Perform following actions correspond to the situation.<br><br>[Case 1] Marker is placed outside the window because of heartbeat.<br>Change FOV size and bed position to fit the marker within the display area.<br><br>[Case 2] Fast heartbeat.<br>The faster the marker moves by heartbeat, the more the marker detection rate goes down.<br>Perform StentShot radiography with higher exposure rate.<br><br>[Case 3] Large FOV size.<br>The larger the FOV size gets, the smaller the marker size gets and detection rate goes down.<br>Set FOV size smaller or set ROI within the range of movement of the marker.<br><br>[Case 4] Large C-arm angle.<br>The deeper the C-arm angle gets, the more the marker contrast and detection rate goes down.<br>Reach performance limits. Perform StentShot radiography with shallower angle. |

| Phenomenon   |
|--|
| When showing StentShot image, delay a showing of next/previous image.<br>Slow reaction and operation of a mouse.<br>Image does not transfer from ACQ to REF. |
| Assumed Cause  |
| PC is highly-load excessively.   |
| Action   |
| Do not perform StentShot radiography for a while (a minute or two).<br>If it does not improve for a long time, restart the system.                           |

|  |
|--|
| Phenomenon   |
| StentShot does not activate, even when pressing a hand switch.   |
| Assumed Cause  |
| Software error, hardware malfunction.  |
| Action   |
| Restart the system. If it is not improved after restarting the system, perform normal fluoroscopy and radiography, and contact our service agency. |

|  |
|--|
| Phenomenon   |
| Abnormal image of StentShot radiography is shown.  |
| Assumed Cause  |
| StentShot parameter is not setting correctly.  |
| Action   |
| Select other DUP once and then select StentShot DUP again.<br>If abnormal image is shown for normal fluoroscopy, restart the system. |

|   |
|---|
| Phenomenon  |
| Result does not display correctly on the StentShot area though setting ROI. |
| Assumed Cause   |
| Marker is placed outside the ROI.   |
| Action  |
| Set ROI including the range of movement of marker by heartbeat.             |

This page is intentionally left blank.

# *Chapter 13*

---

## SCORE Chase (Option)

### ■ Description

|                            |      |
|----------------------------|------|
| 13.1 Overview .....        | 13-2 |
| 13.2 Components .....      | 13-3 |
| 13.3 Specification .....   | 13-5 |
| 13.4 Handling Method ..... | 13-7 |

# 13.1 Overview

## 13.1.1 Purpose

---

SCORE Chase Option (hereafter called "SCORE Chase") is a feature attached to Shimadzu Digital Radiography System DAR-9500f. Perform radiography by moving the table to combine acquired images of different parts of the body, and create a stitched image.

SCORE Chase is useful when treat over a wide range of body parts such as the lower extremity.

## 13.1.2 Feature

---

SCORE Chase has the following features.

- Perform radiography by moving the table, cut and combine small areas filled with contrast agent, and generate and display a stitched image.
- When combined with KS-100, move the table remotely and enable to generate a stitched image without X-ray operator to be exposed to radiation. And perform radiography twice with the same table movement to generate a stitched DSA image which extracted only a blood-vessel image.
- When combined with KS-100, enable to move the table to selected position on the stitched image to reduce X-ray exposure.
- When combined with KS-100, perform an additional radiography at specified position to create DSA image.
- Use an acquired image to change frames for use and reference height, and generate a stitched image with post-processing.
- Display current X-ray exposure position on the stitched image and set radiation position to desired treatment region without being exposed to radiation.



## 13.2 Components

For the system with Cardiovascular Catheterization Table KS-100, the Chase Console is attached to SCORE Chase option.


The table can operate remotely from the control room by using the Chase Console.



System Display

| No. | Part Name     | Quantity |
|-----|---------------|----------|
| 1   | Chase Console | 1 set    |

### 13.2.1 Chase Console



|                 |  |
|-----------------|--|
| Function Button | Use the function button when assigned to enter the RemoteChase mode and move the table to initial position.<br>Table position at the time of pressing the key is registered as radiography start position of RemoteChase and RemoteChase-DSA.<br>And after the test run mode and RemoteChase-DSA Live acquisition, continue to press the key to move the table to the start position of radiography. |
|-----------------|--|

### 13 SCORE Chase (Option)

---

|               |  |
|---------------|--|
| Speed Trigger | Use this trigger for table operation during RemoteChase and test run.<br>The table speed will change automatically depending on how hard you pull the trigger. |
|---------------|--|

# 13.3 Specification

## 13.3.1 Specification of SCORE Chase

|                  |  |
|------------------|--|
| Acquisition Rate | <p>&lt;ManualChase /RemoteChase&gt;</p> <p>[DA]</p> <p>Single Radiography: 30 fps, 15 fps, 10 fps, 7.5 fps</p> <p>Bi-plane Radiography: 15 fps, 10 fps, 7.5 fps</p> <p>[RSM]</p> <p>SFD-1212</p> <p>Single Radiography: 12 fps, 6 fps</p> <p>Bi-plane Radiography: 6 fps</p> <p>SFD-1612AF/SFD-1212AF</p> <p>Single Radiography: 15 fps, 7.5 fps</p> <p>Bi-plane Radiography: 7.5 fps</p> <p>SFD-0808/SFD-0808AF</p> <p>Single Radiography: 15 fps, 10 fps, 7.5 fps</p> <p>Bi-plane Radiography: 15 fps, 10 fps, 7.5 fps</p> <p>&lt;RemoteChase-DSA/SPOT-DSA&gt;</p> <p>SFD-1212</p> <p>Single Radiography: 6 fps</p> <p>Bi-plane Radiography: 6 fps</p> <p>SFD-1612AF/SFD-1212AF/ SFD-0808/SFD-0808AF</p> <p>Single Radiography: 7.5 fps</p> <p>Bi-plane Radiography: 7.5 fps</p> <p>&lt;SPOT-DSA&gt;</p> <p>SFD-1212</p> <p>Single Radiography: 6 fps</p> <p>SFD-1612AF/SFD-1212AF/ SFD-0808/SFD-0808AF</p> <p>Single Radiography: 7.5 fps</p> |
| Image Size       | <p>[Acquisition Image]</p> <p>1024 x 1024 Matrix</p> <p>[Stitched Image]</p> <p>Max. 4096 x 16384 Matrix</p>   |

### 13.3.2 Features

#### 13.3.2.1 Image Acquisition

| No. | Items           | Features   |
|-----|-----------------|--|
| 1   | ManualChase     | Perform radiography by moving the table freely without pre-registered table position information, and combine acquired images to generate a stitched image.  |
| 2   | RemoteChase     | Perform radiography by moving the table remotely with the Chase Console Speed Trigger, and combine acquired images to generate a stitched image. As remote operation becomes available, enable to reduce X-ray exposure to physicians.   |
| 3   | RemoteChase-DSA | Inject contrast agent for the first radiography and acquire an image with RemoteChase. And then perform second radiography (mask acquisition) on the same trace to generate a stitched DSA image which extracted only a blood-vessel image.  |
| 4   | SPOT-DSA        | Specify the image acquired with ManualChase or RemoteChase to position the C-arm and table, and perform additional radiography to create DSA image. And then use the last frame of acquired image to enable DSA-MAP automatically. As creating DSA image with acquired image, enable to reduce contrast agent. |


#### 13.3.2.2 Image Processing

| No. | Items              | Features  |
|-----|--------------------|---|
| 1   | Positioning        | According to current table position and height, C-arm position, SID, and FOV, display X-ray radiation position on a stitched image. |
| 2   | PostProcess Stitch | Use an image acquired with Chase mode, change the number of frames used and reference height to regenerate a stitched image.        |

# 13.4 Handling Method

## 13.4.1 Introduction

Select [Chase] or [Chase Console] button, acquisition mode changes to the Chase mode for the next radiography. The Chase button is valid only when the program with enabled Chase mode is selected on Digital User Program (hereafter called "DUP").



## 13.4.2 Manual Chase

Perform radiography by moving the table freely without pre-registered table position information, and combine acquired images to generate a stitched image.



FOV size should be more than 8 inches for radiography.



When operating the C-arm and table tilt or roll during radiography, the image acquired after the operation is not included in a stitched image.




Move the table that contrast agent density peak to be set in the center of image. If not, it may not generate stitched image properly.



Generate a stitched image based on the latest image acquired on each table position. If return the table to position where no contrast agent is filled, blood vessel image may be ended at that position.


- 1 Select a radiography program which enables the Chase mode from radiography program.
- 2 Select [Chase] button on the side menu or SMART Touch.



Side Menu




IVR NEO



SMART Touch



Unable to click [Chase] button while selecting a program which Chase feature is not assigned.



- 3 Modify the reference height from the side menu if needed.  
Select from [ISO Center], [Table Height], [Table Height +5 cm], [Table Height +10 cm] and [Table Height +15 cm].

- 4 Completely filled with contrast agent, and then start radiography.  
Move the table that contrast agent density peak to be set in the center of the image.



If not filled with contrast agent completely and start radiography, contrast blood vessel may not be included in the stitching image. The time from injection of contrast agent to start of radiography is adjustable.

- 5 After radiography, a stitched image is displayed on the reference monitor.



Enable to use measurement and QCA functions for stitched image, but they cause several errors. Therefore, do not diagnose or treat based solely upon this analysis result. Use the analysis result as a reference value.

### 13.4.3 RemoteChase



RemoteChase is enabled when cardiovascular catheterization table KS-100 is installed.



Before RemoteChase, move the table over the entire range for the safety.



When extending the cable for Chase console operation, pay attention not to put pressure on the base of cable unit.



To start RemoteChase after recovering from emergency stop, check if the table moves horizontally by manual or automatic for system status check. If not checking, a message of [Interlock!] is displayed and cannot start RemoteChase.

Perform radiography by moving the table remotely with the Chase console speed trigger, and combine acquired images to generate a stitched image. As remote operation becomes available, enable to reduce X-ray exposure to physicians.

- 1 Select a radiography program which enables the Chase mode from radiography program.
- 2 Move the table to the start position of radiography.
- 3 Press the Chase console function button or [Chase] switch of function switch on the local console.  
[Chase] button on the side menu becomes ON mode automatically.



When assigned to a function switch on the local console, please contact our service representative.



If selecting a program which is not assigned a Chase mode function, cannot press [Chase] button with Chase console function button or function switch on the local console.

- 4 A message of "Start acquisition" is displayed on the acquisition monitor.



Pull the Chase console speed trigger to perform test run and check the table position.

Release the speed trigger to finish test run.

After test run, press the Chase console function button or press [SET] switch of local console or remote console to move the table to the start position.



**NOTE**

To reposition after test run, press [CANCEL] switch to cancel the Chase mode once, and then reposition.

**CAUTION**

When operating the table with Chase console, operate where you can reach to the [STOP] switch of remote console, CyberConsole or CyberGrip immediately.

If the table does not stop though releasing the speed trigger, press [STOP] switch immediately. The table stops.

- 5 Change the reference height of generating a stitched image if necessary.
- 6 Inject contrast agent and start radiography.
- 7 Pull the speed trigger to move the table.

**CAUTION**

When operating the table with Chase console, operate where you can reach to the [STOP] switch of remote console, CyberConsole or CyberGrip immediately.

If the table does not stop though releasing the speed trigger, press [STOP] switch immediately. The table stops.

**NOTE**

If the table exceeds to contrast agent, press [SET] switch on the remote console to move the table in the opposite direction.

- 8 After radiography, a stitched image is displayed on the reference monitor.

### 13.4.4 RemoteChase-DSA

---



RemoteChase-DSA is enabled when cardiovascular catheterization table KS-100 is installed.



Before RemoteChase-DSA, move the table over the entire range for the safety.



When extending the cable for Chase console operation, pay attention not to put pressure on the base of cable unit.




To start RemoteChase-DSA after recovering from emergency stop, check if the table moves horizontally by manual or automatic for system status check. If not checking, a message of [Interlock!] is displayed and cannot start RemoteChase-DSA.

RemoteChase-DSA is an option to generate a stitched DSA image. Inject contrast agent for the first radiography and acquire an image with RemoteChase. And then perform second radiography (mask acquisition) on the same trace to generate a stitched DSA image which extracted only a blood-vessel image.

- 1** Select a radiography program which enables the RemoteChase-DSA from radiography program.
- 2** A message of “Press Function button on Chase Console” is displayed on the top of acquisition monitor.
- 3** Move the table to the start position of radiography.

- 4 Press Chase console function button or [Chase] switch of function switch on the local console.



- 5 A message of “Start acquisition” is displayed on the acquisition monitor.



Pull the Chase console speed trigger to perform test run and check the table position.

Release the speed trigger to finish test run.

After test run, press the Chase console function button or press [SET] switch of local console or remote console to move the table to the start position.



To reposition after test run, press [CANCEL] switch to cancel the Chase mode once, and then reposition.



To start radiography with appropriate X-ray conditions, perform fluoroscopy approximately 2 seconds before radiography.

**CAUTION**

When operating the table with Chase console, operate where you can reach to the [STOP] switch of remote console, CyberConsole or CyberGrip immediately.

If the table does not stop though releasing the speed trigger, press [STOP] switch immediately. The table stops.

- 6 Change the reference height of generating a stitched image if necessary.
- 7 Inject contrast agent and start radiography.
- 8 Pull the speed trigger to move the table.

**CAUTION**

When operating the table with Chase console, operate where you can reach to the [STOP] switch of remote console, CyberConsole or CyberGrip immediately.

If the table does not stop though releasing the speed trigger, press [STOP] switch immediately. The table stops.

**NOTE**

If the table exceeds to contrast agent, press [SET] switch on the remote console to move the table in the opposite direction.

- 9 After radiography, a stitched image is displayed on the reference monitor.
- 10 A message of “Press Function Button to set table to start position” is displayed on the acquisition monitor.  
Press the Chase console function button or press [SET] switch of local console or remote console to move the table to the start position.
- 11 After the table moved to the start position, a message of “Start acquisition” is displayed on the acquisition monitor.  
Press the radiography switch to start mask acquisition (the table moves automatically).
- 12 When the table moves to the same position as live acquisition, automatically complete the acquisition.
- 13 After the acquisition, stitched DSA image is displayed on the reference monitor.

### 13.4.5 SPOT-DSA



SPOT-DSA can be used when Cardiovascular Catheterization Table KS-100 is installed.



SPOT-DSA is available only for frontal radiography.


On SPOT-DSA, specify the image acquired with ManualChase or RemoteChase to position the C-arm and table, and perform additional radiography to create DSA image. And then use the last frame of acquired image to enable DSA-MAP automatically. As creating DSA image with acquired image, enable to reduce contrast agent.

- 1 From Image Selector on the reference monitor, select an image acquired with ManualChase or RemoteChase, and then select a frame of position to be acquired as DSA image.




Cannot apply SPOT-DSA with a different C-arm position and swivel angle, table tilt/roll/swivel angle and FPD rotation angle of image.

- 2 Click [Funct.]-[SPOT] on the side menu of reference monitor and register the target image.



The same DUP (Digital User Program) and FOV as same as the selected image will be selected automatically, and radiography time is changed to 1 second automatically.



- 3 To move to the target position, the C-arm and table positions are notified and the [SET] button blinks. A message of "Press SET key to move to specified position" is displayed on the acquisition monitor. Press each [SET] button and move the C-arm and table to the radiography position.
- 4 After the C-arm and table move to their setting position, a message of "Start acquisition" is displayed. Press the radiography switch and start the mask acquisition.
- 5 After radiography, the last frame is displayed and DSA-MAP mode will be enabled.



If PCIView is enabled, DSA-MAP mode cannot be enabled.

### 13.4.6 PostProcess Stitch


Use the image acquired with Chase mode and regenerate a stitched image.

It is possible to recreate a stitched image from an image not acquired with the Chase mode.




Only DA or RSM images with a FOV size of 8 inch or more and a protocol of 6 fps or higher can be recreated from images not acquired with the Chase mode.

- 1 Select an image (acquired with DA or RSM of at least FOV 8 inch and at least 6 fps) that can be recreated a stitched image from an image selector.
- 2 Click [Funct.]-[Restitch] on the side menu.  
[PostProcess Stitch] window is displayed.



- 3 Register the frame range to generate a stitched image.  
Display the start frame and click [Set Frame]. And display the end frame and click [End Frame]. Or directly input the number of [Start Frame] and [End Frame].
- 4 Register the reference height.  
Check the checkbox when using ISO center. Input of value from 0 to 30 is enabled when specifying the height.
- 5 Click [Preview] to display regenerated stitched image.  
[Create] button is enabled, so click [Create] to save the image.

- 6 Click [Create] button to save the image displayed on a preview as a stitched image, and add to the image selector.



To add an annotation, generate a stitched image first, add an annotation, and then save the image.



### 13.4.7 Positioning

According to current table position and height, C-arm position, SID, and FOV, display X-ray irradiation position on a stitched image. When irradiation position changes due to the operation of table top and C-arm, irradiation position on the stitched image also follow to move.



Not following the C-arm swivel and rotation angle, table tilt, roll and swivel angle, and collimator position.




Display of radiation position is just a simulation result, so it contains error with actual radiation position.

Perform fluoroscopy before acquisition, and then check radiation position.

- 1 Select a stitched image from image selector.
- 2 Click [Funct.]-[Positioning] on the side menu.



For a setting procedure of [Position] button, see ["13.4.10 Setting of SCORE Chase" P.13-23](#) for details.



When moving the table or radiation position disappeared from the window due to zoom-in of the image, an arrow is displayed in the direction of radiation position.






Radiation position is hidden automatically when swivel the C-arm. Radiation position for the stitched image acquired with different C-arm swivel position.

### 13.4.8 Change of Display

On ACQ and REF monitor, enable to change the display to the original image by canceling DSA processing applied with RemoteChase-DSA or SPOT-DSA.

Click [Change] button on the side menu to change the image.


\*If [Change] is assigned to keyboard shortcut and IVR NEO/SMART Touch button, enable to change the display with a shortcut key.



Side Menu




IVR NEO




SMART Touch

### 13.4.9 Image Loop and Identification of Stitched Image

Both image loop and stitched image are displayed as an image selector icon on the image selector. A small icon is displayed on the upper-right of each image selector icon to show whether it is a loop or a stitched image.

| No. | Icon  | Description  |
|-----|---|--|
| 1   |  | An image generated as a stitched image. "S" which indicates as a loop number and a stitched image is displayed on the image selector icon. |

And when a single local study is selected in the study list in the Studies Management window, all series for the study are shown in the Study Series with one icon per type. Series type icons include a brief label beneath them, indicating the type preceded by the number of items within that series.



Possible series type and its corresponding icon is as follows:

| No. | Icon  | Series Type    |
|-----|---|----------------|
| 1   |  | Stitched image |

13.4.10Setting of SCORE Chase

■ Default Setting of Table Position Display

Configuration

System

DICOM Hosts

DICOM

Software

**Hardware**

Hardware (Shutter)

Storage (Local)

Storage (Network)

Display

Menus

Notifications

Physicians

Study Information

Fonts & Colors

Database

Devices

External Software

Fusion

Logs

Hardware Configuration

Information

Grid

FPD Size (inches)

FOD (cm)

ISO Center To Skin (cm)

Dosemeter Configuration

Dosemeter

Table Configuration

☒ Table Reposition

☐ Show ROI indicating the table position

Stroke Length (mm)

Large Monitor Configuration

☒ Large Monitor Control

LMM IP

Field Of View (inches)

|          |                                 |         |                                  |
|----------|---------------------------------|---------|----------------------------------|
| Large    | <input type="text" value="12"/> | FOV 5th | <input type="text" value="4.5"/> |
| Medium   | <input type="text" value="10"/> | FOV 6th | <input type="text" value="0"/>   |
| Small    | <input type="text" value="8"/>  | FOV 7th | <input type="text" value="0"/>   |
| Smallest | <input type="text" value="6"/>  | FOV 8th | <input type="text" value="0"/>   |

Curves Acquisition

☒ Channel 1

☐ Channel 2

☐ Channel 3

☐ Channel 4

☒ Default Acquire Waves

|          |                                  |
|----------|----------------------------------|
| Channel  | Label                            |
| ECG      | <input type="text" value="ECG"/> |
| Pressure | <input type="text" value="ADP"/> |
| PHYSIO   | <input type="text"/>             |
| FLOW     | <input type="text"/>             |

SMART Touch Configuration

Unit

☒ Touch 1

☒ Touch 2

☒ Touch 3


OK

Cancel

Apply

| No. | Item                                   | Description  |
|-----|--|--|
| 1   | Show ROI indicating the table position | When put a check mark, Positioning mode becomes active as a default. |

## Setting IVR NEO and Shortcut Keys




SCORE Chase setting of Shortcut keys tab is configured as follows.




| Category | Button Name | Function  | Icon | A                     | R                     |
|----------|-------------|---|------|-----------------------|-----------------------|
| Actions  | Chase       | Change Enable/Disable of Chase mode.<br>ON= Chase is enabled. |      | <input type="radio"/> | <input type="radio"/> |

A: relevant to Acquisition monitor, R: relevant to Reference monitor

Setting of SMART Touch



SCORE Chase setting of SMART Touch tab is configured as follows.

| Category  | Button Name | Function  | Icon  | A                     | R                     |
|-----------|-------------|---|---|-----------------------|-----------------------|
| Operation | Chase       | Chase mode is available.                          |  | <input type="radio"/> | <input type="radio"/> |
|           | Positioning | Display radiation position on the stitched image. |  | <input type="radio"/> | <input type="radio"/> |
|           | SPOT-DSA    | SPOT-DSA mode is available.                       |  | <input type="radio"/> | <input type="radio"/> |

A: relevant to Acquisition monitor, R: relevant to Reference monitor

This page is intentionally left blank.



# Chapter 14

---

## SMART Display

This chapter describes how to use the Mini Viewer to display images and loops on other computers.

### Description

|      |                                |       |
|------|--------------------------------|-------|
| 14.1 | Introduction .....             | 14-2  |
| 14.2 | Handling Method .....          | 14-4  |
| 14.3 | Setting of SMART Display ..... | 14-12 |
| 14.4 | Trouble Shooting .....         | 14-13 |


# 14.1 Introduction




Available only when SMART Display (option) is installed.

Set SMART Display (option) from reference monitor. Display the monitor setting window according to the following instructions.

- 1 Select [Funct.]-[Monitor] from the side menu.



- 2 Monitor selection window is displayed.





When display [Large Monitor Configuration] window and error occurs, click [Funct.]-[Monitor] again.

If it does not work, please contact our service representative.




## 14.2 Handling Method


### 14.2.1 Main Screen

Enable to select preset, switch video input and customize preset with the mouse and keyboard.



Move the mouse pointer to the right edge of screen to display [Base settings].

Click  (setting language) button to change the language.



#### (1) Base settings

Set language and restart video manager.

#### (2) Preset


Display the preset list displaying on 58-inch monitor. Select and customize preset.

#### (3) Video Input





Display video input displaying on 58-inch monitor.

14.2.2 Basic Settings

Move the mouse pointer to the right edge of screen to display [Base settings].



The following table shows description of each button.

|   |   |
|---|---|
|    | Display general information about the system.                   |
|    | Language of the dialog window/menu will be changed.             |
|    | The dialog window for rebooting SMART Display is displayed.     |
|  | The dialog window for shutting down SMART Display is displayed. |



NOTE

Do not shut down SMART Display.  
Or the power supply of digital acquisition cabinet should be turned on.




NOTE

It takes about 3 minutes to restart SMART Display. Image will not display and operate on 58-inch monitor while restarting.









### 14.2.3 Preset

Move the mouse pointer to the left edge of screen to display [Preset]. Enable to select preset and customize preset.



Left-click the mouse to select preset. The color of currently selected preset is in reverse.

The following table shows description of each button.

|   |  |
|---|--|
|  | Create new preset.   |
|  | Change the name of preset.   |
|  | Set preset a default. Restart 58 inch monitor to display the default preset. |
|  | Edit preset.   |
|  | Duplicate preset.  |
|  | Delete preset.   |
|  | Save the change of preset.   |
|  | Reverse the operation before saving the changes.                             |

14.2.4 Video Input


Move the mouse pointer to the upper edge of screen to display [Video Input].




Thumb nail of video input displaying on 58-inch monitor is displayed.

14.2.5 Switch Segment

- 1 Select the segment to change displaying position.



- 2 Drag and drop to the segment which to be displayed. Segment will be switched.



## 14.2.6 Switch Video Input


- 1 Move the mouse pointer to the upper edge of screen to display [Video Input].



- 2 Select video input to be displayed.



- 3 Drag and drop to displaying segment.





### 14.2.7 Creating Preset


Create new preset according to the following instructions.

- 1 Move the mouse pointer to the left edge of screen to display [Preset].




Select  button.

- 2 [Create new layout] is displayed.




Enter an optional name for the name and click [OK].

- 3 Screen will be displayed in 1 screen (segment).



Click  button.

- 4 [Split window] is displayed. Select desired layout.



- 5 Assign an optional video input to each segment. (Refer to "14.2.6 Switch Video Input" P.14-8)



- 6 After preset is completed, click  button in preset. Preset will be saved.



# 14.3 Setting of SMART Display

## Setting of IP Address

**Configuration**

**System**

DICOM Hosts

DICOM

Software

**Hardware**

Hardware (Shutter)

Storage (Local)

Storage (Network)

Display

Menus

Notifications

Physicians

Study Information

Fonts & Colors

Database

Devices

External Software

Fusion

Logs

**Hardware Configuration**

**Information**

Grid: IN

FPD Size (inches): 12

FOD (cm): 72

ISO Center To Skin (cm): 15

**Field Of View (Inches)**

|          |    |         |     |
|----------|----|---------|-----|
| Large    | 12 | FOV 5th | 4.5 |
| Medium   | 10 | FOV 6th | 0   |
| Small    | 8  | FOV 7th | 0   |
| Smallest | 6  | FOV 8th | 0   |

**Dosemeter Configuration**

Dosemeter: VacuDAP

**Table Configuration**

☒ Table Reposition

☐ Show ROI indicating the table position

Stroke Length (mm): 1350

**Large Monitor Configuration**

☒ Large Monitor Control

LMM IP: 192.168.100.37

**Curves Acquisition**

☒ Channel 1

☐ Channel 2

☐ Channel 3

☐ Channel 4

☒ Default Acquire Waves

| Channel  | Label |
|----------|-------|
| ECG      | ECG   |
| Pressure | ABP   |
| PHYSIO   |       |
| FLOW     |       |

**SMART Touch Configuration**

☒ Touch 1

☒ Touch 2

☒ Touch 3

| Unit        |
|-------------|
| TOUCH-01-01 |
| TOUCH-01-02 |
| TOUCH-01-03 |


OK Cancel Apply

| Item   | Description                      |
|--------|----------------------------------|
| LMM IP | Set IP address of SMART Display. |

# 14.4 Trouble Shooting

## 14.4.1 Error Message

[Error Message] will be displayed at the right bottom of screen if any failure occurs on some hardware.




Blinking red frame indicates hardware with any trouble occurs.


| Display                      | Explanation                           | Action                              |
|------------------------------|---------------------------------------|-------------------------------------|
| FPGA Temperature             | Error occurs on FPGA temperature.     | Contact our service representative. |
| OnBoard Temperature          | Error occurs on Hardware temperature. |                                     |
| Case Font Fan/Additional Fan | Error occurs on FAN.                  |                                     |

## 14.4.2 Displaying Image

If an image is not displayed on 58 inch monitor, check color of LED.




| LED       | Explanation                 | Action  |
|-----------|-----------------------------|---|
| Yellow    | No video input.             | Contact our service representative.   |
| Red       | Failure of 58 inch monitor. |   |
| Light Out | Failure of 58 inch monitor. | Power of 58 inch monitor is turned off. Check that the power of 58 inch monitor is turned on. The power switch is on the back of monitor. |




If an image is not displayed on 58 inch monitor in spite of confirming the way above, use the monitor cart or Folding Backup Monitors to confirm an image.



Monitor cart (Example)



Folding Backup Monitor (in use)



Folding backup monitor (not in use)



Do not use both of the Cart and Backup Monitor Kit for a purpose without mounting backup monitor.



Use the Monitor Cart or Folding Backup Monitors only when SMART Display breaks down.  
Especially, keep the Folding Backup Monitors in full-fold position when not in use.

This page is intentionally left blank.



# *Chapter 15*

---

## Temporary License

### ■ Description

|      |                                 |      |
|------|---------------------------------|------|
| 15.1 | Introduction .....              | 15-2 |
| 15.2 | Start a Temporary License ..... | 15-3 |
| 15.3 | End a Temporary License .....   | 15-4 |

## 15.1 Introduction

This chapter describes how to use a temporary license for 90 days of optional application. The following applications have a temporary license.


| Name            | Function   |
|-----------------|--|
| SCORE StentView | Extract a marker position from a real-time image and enhance and display the device by digital image processing. Add some frames and display enhanced stent image.<br>This is useful when place 2 stent and confirm the position of second one.  |
| SCORE StentShot | Extract a marker position from a real-time image, and enhance and display the device by digital image processing. Add some frames and display enhanced stent images.<br>This is useful for confirming the shape of placed stent.   |
| SCORE Chase     | Perform radiography by moving the table without pre-registered table position information, and combine acquired images to generate a stitched image.<br><br>When combined with KS-100, enable to move the table to the position selected on the stitched image to reduce X-ray exposure. |



Hear an explanation of application usage from application instructor or service engineer, and then start a temporary license.

## 15.2 Start a Temporary License


- 1 On REF monitor, click  mark at the bottom-right of the monitor.
- 2 Click [License Info].



- 3 Select an application to start a temporary license and click [Start].

| License Information |                                      |                |                 |               |
|---------------------|--------------------------------------|----------------|-----------------|---------------|
|                     | Feature                              | Mode           | Expiration Date | Remaining Day |
| Frontal             | DAR9500                              | Permanent      |                 |               |
| Frontal             | StentView_Temporary_EndUser          | Not Started    | 0/0/0           | 0             |
| Frontal             | StentView_Temporary_Installation     | Limited period | 2018/2/27       | 14            |
| Frontal             | StentShot                            | Permanent      |                 |               |
| Frontal             | StentShot_Temporary_EndUser          | Not Started    | 0/0/0           | 0             |
| Frontal             | AutoStitching_Temporary_EndUser      | Permanent      |                 |               |
| Frontal             | AutoStitching_Temporary_Installation | Limited period | 2018/3/14       | 29            |
| Frontal             | FlexAPS                              | Permanent      |                 |               |

- 4 A message dialog of "Do you want to start the temporary application?" is displayed. Click [Yes] to start a temporary license. Application is available for 90 days.



Expiration date and the remaining days are displayed every time start the system.

## 15.3 End a Temporary License

An application will end after 90 days from the start of use.

End the temporary license with the following procedure.




**NOTE**

Unable to restart a temporary license once you have completed.

- 1 Select an application to end a temporary license on [License Information], and then click [End].

| License Information |                                      |                |                 |               |
|---------------------|--------------------------------------|----------------|-----------------|---------------|
|                     | Feature                              | Mode           | Expiration Date | Remaining Day |
| Frontal             | DAR9500                              | Permanent      |                 |               |
| Frontal             | StentView_Temporary_EndUser          | Not Started    | 0/0/0           | 0             |
| Frontal             | StentView_Temporary_Installation     | Limited period | 2018/2/27       | 14            |
| Frontal             | StentShot                            | Permanent      |                 |               |
| Frontal             | StentShot_Temporary_EndUser          | Not Started    | 0/0/0           | 0             |
| Frontal             | AutoStitching_Temporary_EndUser      | Permanent      |                 |               |
| Frontal             | AutoStitching_Temporary_Installation | Limited period | 2018/3/14       | 29            |
| Frontal             | FlexAPS                              | Permanent      |                 |               |

- 2 A message dialog of "Do you want to end the temporary application?" is displayed.



Click [Yes] to end.

# *Chapter 16*

---

## Mini Viewer


This chapter describes how to use the Mini Viewer to display images and loops on other computers.

### Description

16.1 The Mini Viewer .....16-2








## 16.1 The Mini Viewer







When writing a DICOM CD, the Mini Viewer application (EIVIEWER.exe) is included with the images files on the DICOM CD. The Mini Viewer can be used to view loops and images on any Windows computer. When you launch the Mini Viewer (double-click EIVIEWER.exe), the first image loops is displayed like this:



The following information is displayed at the edges of the window:

| Location | Description   |
|----------|---|
| Top      | Patient name and study information.                                   |
| Right    | Image Selector, with image icons.                                     |
| Bottom   | Playback / Pause button and Brightness / Contrast adjustment buttons. |
| Left     | Function buttons.   |

| No. | Button/GUI  | Name               | Function  |
|-----|---|--------------------|---|
| 1   |  | Playback/<br>Pause | Click once to start playback and again to pause.  |
| 2   |  | Study              | Display the list of studies from which you can select one to review.  |
| 3   |  | Zoom               | Click [Zoom] once to make the complete image fit into the image area. Click [Zoom] again to display the image to actual size.   |
| 4   |  | BMP                | Save the displayed image (pause loop first) in a Windows BMP (bitmap) image file.   |
| 5   |  | Info               | Display version information (About tab) and Help information (Help tab).  |
| 6   |  | RDSR               | Click to check the description of RDSR.<br>See  <a href="#">"6.2.4 Displaying RDSR" P.6-13</a> |

| No. | Button/GUI  | Name                | Function   |
|-----|---|---------------------|--|
| 7   |    | Brightness/Contrast | With the mouse, drag the slider (vertical line) to the right to increase Brightness / Contrast, or to the left to decrease it.                                       |
| 8   |    | Full-screen Mode    | Click the small square in the upper-right window corner to display switch to full-screen mode. Left-click anywhere in the Mini Viewer to switch back to normal size. |
| 9   |    | Window size         | You can make the Mini Viewer window larger by dragging the size handle down and to the right, or reduce its size by dragging the handle upward and to the left.      |
| 10  |    | L n.n               | Series number. Loop number   |
|     |   | nnn/nnn             | Frame number / total number of frames  |
|     |   | n.nnx               | Zoom level, 1.00 =100 %, 1.50 = 150 % of original size.  |
|     |   | nnF/S or Pause      | Playback speed in frames per second, or "Pause".   |
| 11  |   | Image Selector      | Click the icon of the desired loop or Scroll the Image Selector up or down by right-clicking and dragging the Image Selector.  |
| 12  |  | Exit Mini Viewer    | Click the "X" in the upper-right corner to close the Mini Viewer.  |



#### Keyboard Operation

The keyboard arrow keys can be used to control playback. Press Left / Right to move to the previous / next loop image within the selected loop. Press Up / Down to move to the previous / next loop. Press the space bar to switch between playback and pause.



Biomedical signals such as ECG curves and parameters such as FOV and SID, are not displayed by the Mini Viewer.

### ■ Install DirectX 11

The Mini Viewer requires at least DirectX 11.



Do not attempt to install DirectX 11 or use the Mini Viewer on a DAR-9500f computer. It is intended only for use on other Windows computers.




Consult Microsoft documentation for full installation instructions. The information provided here is in summary form.

- **Check Current Version**

First if at least DirectX11 is already installed.


Set up according to the following procedure (for Windows7).

- 1 On the Windows Start menu, choose [Run].
- 2 Enter "dxdiag" and click [OK].





The DirectX Diagnostic Tool appears. Check the DirectX version item near the bottom (marked with a dashed box in the image below).



If this dialog does not appear, DirectX might not be installed or an old version might be present. Refer to the following paragraph to get DirectX 11 and to install or update it.

- **Installing DirectX 11**

If you do not have at least DirectX 11, follow the Microsoft instructions on how to get and install it at:

<http://www.microsoft.com/ja-jp/download/details.aspx?id=35>



## CAUTION

Please note that we assume no responsibility for any problems encountered as a result of installing DirectX 11.

This page is intentionally left blank.

# Chapter 17

## Administration

This chapter describes how to administer the system and includes the following.

### Description

|      |                                   |       |
|------|-----------------------------------|-------|
| 17.1 | Introduction .....                | 17-2  |
| 17.2 | Switching Modes .....             | 17-2  |
| 17.3 | Viewing User Type .....           | 17-3  |
| 17.4 | Storage Commit Messages .....     | 17-4  |
| 17.5 | Working with LUT Profiles .....   | 17-4  |
| 17.6 | DICOM Monitor Linearization ..... | 17-10 |
| 17.7 | Options Configuration .....       | 17-14 |

## 17.1 Introduction

Administration mode enables you to fully configure the product and access individual network servers and CD/DVD drives.



## 17.2 Switching Modes


You switch between User and Administration mode by logging in to your operating system with a user name that is configured in DAR-9500f to have the desired user type, either user or administrator.

### 17.2.1 Switch to Administration Mode

---

Follow this procedure to switch from User mode to Administration mode:

- 1 Execute the shutdown procedure as described in  ["3.2 User Authentication" P.3-4](#).
- 2 Execute the startup procedure as described in  ["3.1 Startup" P.3-2](#).  
After Windows startup, hold down the keyboard [Ctrl]+[Alt]+[Delete] key to force display of the Windows login prompt and log off.
- 3 Log in with an administrator user name "super".  
Double-click the DAR-9500f icon on the desktop to start the software. The GUI starts in Administrator mode.





Now, when you click the [Options] button in the image Viewer window, you will be able to configure any option.

### 17.2.2 Switch to User Mode

---

To switch back to User mode, follow this procedure:


- 1 Execute the shutdown procedure as described in  ["3.2 User Authentication" P.3-4](#).
- 2 Execute the startup procedure as described in  ["3.1 Startup" P.3-2](#).  
Windows automatically logs in as user and the GUI starts in User mode.

## 17.3 Viewing User Type

Regular users of DAR-9500f should log in to the operating system with a user name that has been configured in DAR-9500f to have a user type of [USER]. Normally, automatic login to User mode is performed. Only the system administrator should log in with a DAR-9500f user type of [SUPER].

Follow this procedure to view and optionally set user types:

- 1 Run DAR-9500f in Administration mode as described above in ["17.2.1 Switch to Administration Mode" P.17-2](#).
- 2 Click [Options] on the side menu.  
The tabbed Configuration dialog box appears.
- 3 Click [Database] tab.  
It looks like this.



- 4 In the upper-right corner, click [User Table].  
The [User] dialog box appears like this:

| User   |            |            |           |
|--|------------|------------|-----------|
|  | User Login | User Name  | Type      |
| D  | Admin      | admin      | SUPER     |
|  | operator   | operator   | USER      |
|  | super      | super      | SUPER     |
|  | serviceapp | serviceapp | SUPER     |
|  | shimadzu   | shimadzu   | SUPER     |
|  | adminope   | adminope   | SUPERUSER |
| <div> Add/Modify Delete Update Refresh Close Set Password </div> |            |            |           |



### CAUTION

Do not change the User Type settings. Otherwise, service personnel will not be able to perform system maintenance or the system may not operate properly.

## 17.4 Storage Commit Messages

If the DICOM archive server supports it, Storage Commit can be used to ensure greater integrity in all server updates. Several messages can appear related to Storage Commit as follows.

- The “Server\_Name” does not support Storage Commit.  
In the [AE Hosts] table, uncheck the [Storage Commit] option for this server.
- The maximum number of retries without answer on server “Server\_Name” for the study “Study\_Name” is exceeded.  
This can occur when the network is disconnected for a long period of time. Normal operation is resumed soon after the network connection is re-established.
- The maximum number of failures is reached on server “Server\_Name.”  
This indicates that there may be a problem with the network connection or file format. Contact the DICOM archive server administrator.

## 17.5 Working with LUT Profiles

### 17.5.1 Introduction

---

Images are acquired with 12 bits (4096 levels) of grayscale, manipulated as 32-bit data, and then displayed or printed as 8-bit data (256 levels of grayscale).

The many gray levels in the image data must be mapped to fewer gray levels for display and printing. By default, a linear mapping scheme is used. However, since most display monitors and printers do not display or print gray levels in a linear manner, custom display and printer profiles can be created to match the non-linear device characteristics. Furthermore, due to differences in lighting, personal preference, and perception, display monitors sometimes need additional adjustment beyond basic linearization.


In DAR-9500f, mapping profiles are defined independently for monitor and printer via what is called a lookup table (LUT).

Printer linearization and display adjustment profile LUTs are defined graphically via a simple curve editor.

The curve is converted into a lookup table (LUT) that is applied to the image data at display and print time.




#### NOTE

Display monitor linearization profiles are defined at installation time according to the [DICOM part 14 Grayscale Standard Display Function]. For details, see  ["17.6 DICOM Monitor Linearization" P.17-10.](#)


## 17.5.2 LUT Profile Definition

To define either printer linearization or display adjustment LUT profiles, follow this procedure.

- 1 If you are not ready in Administration mode, switch to it now as described in ["17.2.1 Switch to Administration Mode" P.17-2.](#)
- 2 Display a reference image, preferably one featuring all possible gray levels.
- 3 Click [Filters] on the side menu.  
The [LUT] fly-out menu appears.



- 4 Click [LUT-O] button.  
The [LUT Editor] dialog box appears.




The curve portion of the [LUT Editor] dialog box, defines the profile graphically as follows:

- The permanent dashed line from the lower-left corner to the upper-right corner represents the input data with a value of 0 to 255. With a linear profile, the output is the same as the input.
- The curve that you draw defines what the output value will be, also in the range 0 to 255. Points above the default dashed line will cause the output to be brighter (higher number) than the input. Conversely, points below the dashed line will cause the output to be darker (lower number) than the input. Whenever the pointer is in the graph area, both the input and output values are shown for the pointer position. In this example, we see that a gray input value of 80 is darkened to a gray output value of 57.

- This sample graph increases contrast as follows. It causes gray levels below middle gray (center of graph) to be made darker (added points are below the dashed line) and gray levels above middle gray to be made lighter (added points are above the dashed line). Note also how that the pure black (0) and pure white (255) inputs are output unaltered as 0 and 255 respectively, because the curve starts and ends on the default linear center line.

- 5 In the drop-down list to the right of the [Selection] button, choose either [Lut Custom] to define a display profile or [Lut Print] to define a printer profile.




- 6 If an existing curve is similar to the curve you wish to define, click [Load] and choose it as your new curve starting point.

- 7 In the curve portion of the [LUT Editor] dialog box, define your curve by clicking points in the graph area. Connecting lines are drawn automatically from the curve box corners and between the points. Reposition a point by dragging it. Delete a point by clicking it with the left and right mouse buttons pressed simultaneously.

For display profiles, the image appearance changes as you shape the curve. For all profile types, the input and output values are updated as you move the pointer.

- 8 When satisfied with the curve, click [Save], enter a meaningful name (up to 50 characters), and click [Save].

- 9 To delete a profile, click [Delete] and then the profile name.



You may find it convenient to create several profiles and then see how various images display or print with the profiles applied to them.




17.5.3 Assigning Display Adjustment LUT Profiles

Once you have defined a display adjustment LUT profile, you can assign it to one of 10 positions in the profile list as follows.


- 1 In the drop-down list to the right of the [Selection] button, choose [Lut Custom].




- 2 Click [Load] and then the profile that you wish to assign.




- 3 Click [Selection] and then click the position in the list to which you wish to assign this profile.




When in User mode, the assignments that you make here are shown as choices on the [Filter]-[LUT-O] list.



- 4 To delete an assignment, click [Selection] and then the assignment that you wish to delete. Assignments have a check mark to their left.



It is OK to adjust the three predefined profiles LUT-A, LUT-B and LUT-C, however, do not change the shape of these curves so much that they no longer resemble the shapes shown on the buttons.




### 17.5.4 Assigning Printer LUT Profiles

Once you have defined a printer LUT profile, you can assign it to one of 10 positions in the print profile list as follows.


- 1 In the drop-down list to the right of the [Selection] button, choose [LUT Print].



- 2 Click [Load] and then the profile that you wish to assign.




- 3 Click [Selection] and then click the position in the list to which you wish to assign this profile.



The assignments that you make here are shown as choices in the LUT list of the print windows. See the next section ["17.5.5 Setting the Default Printer LUT Profile" P.17-9](#).

- 4 To delete an assignment, click [Selection] and then the assignment that you wish to delete.


Assignments have a check mark to their left.



### 17.5.5 Setting the Default Printer LUT Profile

To set the default printer LUT profile, follow this procedure.

- 1 If you are not already in Administration mode, switch to it now as described in ["17.2.1 Switch to Administration Mode" P.17-2.](#)
- 2 Click [Funcnt] and [Print] on the side menu.



- 3 The print window appears. Click the [LUT] button.  
A small profile selection box appears at the upper-right corner of the print window.



- 4 Choose the desired default profile from the drop-down list.  
The chosen profile name appears near the LUT button. If desired, click [Close] to close the drop-down list.
- 5 To test the selected profile, click [Add] to add the current image to the print list with the selected profile applied.  
You can preview and print the same image with several different profiles applied and then choose the best one as the default. When finished, click [Close] to close the print window.

## 17.6 DICOM Monitor Linearization



This section gives basic information on how to perform DICOM monitor linearization using the tool available on the LUT Editor dialog box. You must provide your own photometer and are encouraged to consult the DICOM PS3.14 standard.

Typically performed by the installer, the Reference monitor can be linearized with a photometer and the provided test pattern. Luminance readings can be taken from the monitor for the full range of the video adapter Discrete Data Levels (DDLs). The readings can then be entered in a grid to compute the characteristic display curve LUT for the monitor.


Luminance is defined as the luminous intensity per unit area projected in a given direction, measured in candela per square meter ( $\text{cd}/\text{m}^2$ ). The luminance generated by an emissive display system may be measured with a photometer.

DDL (Discrete Data Level) is defined as a digital value, typically in the range 0 to 255 (0=darkest, 255=lightest) (for each of the 3 colors red, green, and blue) sent to the video adapter to produce luminance levels on the monitor.

Follow this procedure to perform DICOM monitor linearization. (Recommended photometer: Siemens, SMfit ACT Calibration).

- 1 Log in to the DAR-9500f.
- 2 Load an image or loop that contains as many gray levels as possible.
- 3 Wait for the monitor to warm up for at least 10 minutes.
- 4 Set the room lighting to a normal level.
- 5 Choose [Filters]-[LUT-Others], and then click [DICOM] button in the lower-right corner of the [LUT Editor].

The [LUT Editor] dialog box expands to the right with the DICOM Linearization panel.






Unless otherwise started, all the following button references are for the [DICOM Linearization] panel.

- 6 Use the [Delete] button to delete any existing values in the [Monitor Characteristic Curve] grid.
- 7 Click the [Test Pattern] button.  
The test pattern is displayed with a 20% gray background and a square (defaults to white) in the middle.
- 8 If not already set to 255, adjust the DDL value of the square to 255 (100%) by pressing and holding the right-arrow key.
- 9 Using the Brightness/Contrast controls of the monitor set the Brightness and Contrast so that the white square in the middle is the brightest white possible without making the background gray (20% gray) too light.
- 10 Place the photometer on the center of the monitor, over the white square.
- 11 With the DDL value at 255 (100%), record the luminance value reported by the photometer.
- 12 Decrease the DDL value by pressing the left-arrow key until a change in luminance is detected.
- 13 Record the DDL value and luminance.
- 14 Repeat this process until you have at least three DDL/Luminance pairs. Normally, at least ten pairs are needed for good curve definition.  
The readings should be spaced over the entire 0-255 DDL range.
- 15 Enter the largest DDL value and its corresponding Luminance value in the two boxes below the grid and click [Add]. If, before clicking [Add] you wish to clear both entry boxes, click [Clear].
- 16 Add all other DDL/Luminance pairs in descending-DDL order.




The luminance value entered MUST BE less than or equal to the Luminance value entered for any larger DDL. For example, if you add DDL/Luminance pair 215/305.2890, the Luminance value for DDL 214 must be no more than 305.2890.

- 17** To change a value in an existing row, click the row (an arrow appears at the left edge), make your changes in the two boxes below the grid, and then click [Update].
- 18** To delete a row from the grid, click the row (an arrow appears at the left edge) and then click [Delete].
- 19** Once you are satisfied with your list of DDL/Luminance pairs, click [Compute]. The display characteristic curve is calculated and saved as a LUT.  
A message to this effect is displayed




- 20** If there are problems with your DDL/Luminance pairs, an error message indicating that incoherent values were encountered is displayed.




If this happens, adjust the values so that the rule described in step 16 above is respected and click [Compute] again.

- 21** To apply your new display linearization LUT, click [Apply LUT] so that the check mark appears. The visual effects are immediate.  
This LUT will remain in effect for all future image viewing on this monitor as long as Apply LUT remains checked.



- 22** To close the [DICOM Linearization] panel of the [LUT Editor], click the arrow at the upper-left edge of the [Monitor Characteristic Curve] grid.





If the monitor or video adapter is replaced (even with the same model) redo the DICOM monitor linearization.

## 17.7 Options Configuration

### 17.7.1 Introduction

Although the Installation personnel perform the initial configuration according to your institution's specifications, you can further adjust the configuration as needed. Configuration is performed via the [Options] tab in the [Studies Management] window and the [Options] button on the side menu.

### 17.7.2 Options Tab

Three basic configuration items are available on the Options tab of the Studies Management window. Click the **Options** tab to reveal its three check boxes which are defined as follow.



| No. | Check Box   | Description  |
|-----|---|--|
| 1   | Automatic Refresh (with adjustable minutes value) | When checked (a check mark appears in the check box), the Study List is automatically refreshed (redrawn) on the screen, showing any new studies that have arrived since the last refresh or other status changes. The refresh interval in minutes, is entered in the combo box to the right of [Automatic Refresh]. Permissible values are 1 to 120 minutes with a default of 30.   |
| 2   | "All Time" study query                            | When "All Time" study query is checked, the default Date search criteria automatically entered in the Date box is changed from today's date to [All Time], meaning all studies in the Study List.  |
| 3   | Search by keywords                                | Normally, with [Search by keyword] unchecked, search criteria matching is only performed at the beginning of the searched text so that with the Patient Name criteria set to "rog", only last names <i>beginning</i> with "rog" such as "Rogers" and "Roget" would be found. Names such as "Smith Roger" would not be found. With [Search by keyword] checked, search criteria matching is performed throughout the searched text so that with the same Patient Name criteria of "rog", all patient names with the three letters (in sequence) <i>anywhere</i> in the searched text would be found, including names such as "Smith Roger." |



### 17.7.3 Configuration Dialog Box

To display the Configuration dialog box, click [Options] on the side menu of either the Reference or Acquisition monitors. The tabbed Configuration dialog box appears with the System tab initially displayed.




#### CAUTION

In the option description tables appearing throughout this section, do not modify items with an “X” mark in the “Mod” column. The system might not operate properly if such items are modified. Items marked “OK” can be modified according to supplied information.



#### NOTE

If you only see the Menus, Display, Software, Devices, Notifications, and Fonts and Colors tabs, you are still in User mode. If you are not already in Administration mode, switch to it now as described in  ["17.2.1 Switch to Administration Mode" P.17-2](#).



#### NOTE

If you sign in with user mode, the Cyber Security tab appears only if you sign in at the SUPER level.



#### NOTE

DAR-9500f systems have two computers, one for Reference with its monitor and the other for Acquisition with its monitor. For users, a keyboard and mouse is only available for the Reference computer. To configure the Acquisition computer, the administrator must make sure that a mouse and keyboard is connected to the Acquisition computer (these may already be connected, but hidden from the user). It is not necessary to turn off the Acquisition computer before connecting the mouse and keyboard.



#### NOTE

Unless otherwise directed, enter identical configuration values on both the Reference and Acquisition computers. When indicated, enter values only on the Reference computer. They will be automatically propagated to the Acquisition computer.



In this Manual, when a Configuration dialog box tab is identical for Reference and Acquisition, only one is shown. When there is a difference, both are shown.

System Tab

System options are configured as follows.

Reference Monitor

Configuration

System

DICOM Hosts

DICOM

Software

Hardware

Hardware (Shutter)

Storage (Local)

Storage (Network)

Display

Menus

Notifications

Physicians

Study Information

Fonts & Colors

Database

Devices

External Software

Fusion

Logs

Cyber Security

System Configuration

System Information

Computer Name: DAR9500REV-01

Institution Name: Institution Name

Institution Address: Institution Address

Department Name: Department

Station Name: DAR9500REV-01

Machine Serial Number: 0000000011302

DICOM AE Title: DAR9500REV\_01

Manufacturer Model Name: Trinias

Implementation UID: 1.2.392.200036.9110.18.11302

Information Version: Voyager\_V6\_10\_0

Link

Link Mode: ☒

Linked Station: 192.168.100.22

OK Cancel Apply

Acquisition Monitor

Configuration

System

DICOM Hosts

DICOM

Software

Hardware

Hardware (Shutter)

Storage (Local)

Storage (Network)

Display

Menus

Notifications

Physicians

Study Information

Fonts & Colors

Database

Devices

External Software

Fusion

Logs

System Configuration

System Information

Computer Name: DAR9500ACQ-F-01

Institution Name: Institution Name

Institution Address: Institution Address

Department Name: Department

Station Name: DAR9500REV-01

Machine Serial Number: 0000000011301

DICOM AE Title: DAR9500ACQ\_F\_01

Manufacturer Model Name: Trinias

Implementation UID: 1.2.392.200036.9110.17.11301

Information Version: Voyager\_V6\_10\_0

Link

Link Mode: ☒



Biplane: ☐


(A) Frontal: ☐

(B) Lateral: ☐

Plane B Name:

OK Cancel Apply

| Item                      | Description  | Mod |
|---------------------------|--|-----|
| <b>System Information</b> |  |     |
| Computer Name             | The name of this computer configured in Windows. Enter the name ONLY on the Reference computer. It is propagated to the Acquisition computer.  | ×   |
| Institution Name          | <p>Administrator-defined institution name set according to your needs. Enter the name ONLY on the Reference computer. It is propagated to the Acquisition computer. (Up to 64 characters (no backslash “\”), 35 recommended maximum.) You must shut down and restart the system after changing this item.</p> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;"> <div style="text-align: center;">  <b>NOTE</b> </div> <p>When the system needed activation during installation, system activation will be needed again if the hospital name is changed.</p> <div style="border: 1px solid black; padding: 5px; margin: 10px 0;"> <p><b>Message</b></p> <p> If the hospital name is changed, system activation will be needed again. Do you really want to change the hospital name?</p> <div style="text-align: right;"> <span>Yes</span> <span>No</span> </div> </div> </div> | ○   |
| Institution Address       | Administrator-defined institution address set according to your needs. Enter the address ONLY on the Reference computer. It is propagated to the Acquisition computer.   | ○   |
| Department Name           | Administrator-defined department name set according to your needs. Enter the name ONLY on the Reference computer. It is propagated to the Acquisition computer. (Up to 64 characters (no backslash “\”), but recommended maximum of 35.) You must shut down and restart the system after changing this item.   | ○   |
| Station Name              | Administrator-defined station name. The name entered here appears in the Origin column of the Studies Management window. Enter the name ONLY on the Reference computer. It is propagated to the Acquisition computer. (Up to 16 characters: alpha, number, hyphen “-”, underscore “_”, period “.”, apostrophe “”’.”.)  | ○   |
| Machine Serial Number     | A number (should be unique for this type of system) that identifies this computer and is used in creation of the Study Unique ID. The software must be restarted after changing this option so it automatically exits. (Up to 13-digit number with forced leading zeros if fewer than 13 digits.)  | ×   |
| DICOM AE Title            | The DICOM device name of the computer. The default is auto-generated from the Machine Serial Number. (Up to 16 characters: alpha, number, underscore “_”, period “.”.)   | ×   |
| Manufacturer Model Name   | Enter a product model name (such as Trinias) if necessary.   | ○   |


| Item                | Description  | Mod |
|---------------------|--|-----|
| Implementation UID  | Fixed tags appearing in the DICOM files produced by this product.  | ×   |
| Version Information | Displays Version of this product.  | ×   |
| <b>Link</b>         |  |     |
| Link Mode           | Always enabled, specifies that the Reference and Acquisition computers are linked to form the DAR-9500f system.  | ×   |
| Linked Station      | (Reference only) Must be set to the name shown in the [Computer Name] field on the Acquisition computer Configuration dialog box or IP address of Acquisition computer.  | ×   |
| Biplane             | (Option) (Acquisition only) Enables Biplane mode in which two separate acquisition computers are used for the simultaneous acquisition of images at two angles. Both computers must have Biplane checked. One computer must also have [(A) Frontal] checked and the Windows machine name of the other computer (plane B) entered in [Plane B Name], and the other computer must have [(B) Lateral] checked.<br>See  "10 Difference with Bi-plane Option" P.10-1 | ×   |

## DICOM Hosts Tab

The Hosts tab defines DICOM devices with which this DAR-9500f can communicate. Its options are defined as follows.

| Item                  | Description   | Mod |
|-----------------------|---|-----|
| <b>Information</b>    |   |     |
| AE Title              | The DICOM device name of the station. (Up to 16 characters.) To enter a new host, Click [New] and fill in the information.  | ○   |
| Echo                  | To verify that the connection to a DICOM device such as a server is valid, choose the name in the drop-down list below the Echo button and optionally check [Ping]. Click [Echo] to perform the test. If [Ping] is selected, the device is first pinged, and then the DICOM echo is performed. As appropriate, success or failure messages are displayed. | ○   |
| Host                  | The TCP/IP protocol host name or IP address of the DICOM device. (Up to 64 characters.)   | ○   |
| Alias                 | Alias for Host. When using the Send feature on the Studies Management window, this name appears in the Destination list. Also, when in Administration mode, this name appears in the Location drop-down list. <b>Do not enter more than 11 characters.</b>  | ○   |
| Port                  | The port number used by the DICOM device. Most DICOM devices default to port 104. (0 to 65535.)   | ○   |
| Packet Size           | The maximum packet size supported by the DICOM device. The default of 99280 is usually acceptable. (10,000 to 1,000,000 (one million).)   | ○   |
| File Packets Grouping | Defines how many packets should be grouped together for socket-level communications (1 to 5).   | ○   |
| UID/PSW/DSN           | Unused except for external studies.   | ○   |
| <b>DICOM Host</b>     |   |     |
| Yes/No                | Indicates whether it is a DICOM host or not.  | ○   |

| Item                     | Description   | Mod                      |
|--------------------------|---|--------------------------|
| <b>Host Type</b>         |   |                          |
| Modality Worklist Server | The DICOM device interfaces to a Hospital Information System (HIS). When selected, and a new study is started, the chosen modality worklist server will be queried for procedures waiting to be done.   | <input type="radio"/>    |
| Printer                  | The device is a DICOM networked printer.  | <input type="radio"/>    |
| Other                    | The DICOM device is a server or destination station.  | <input type="radio"/>    |
| <b>Capabilities</b>      |   |                          |
| Store                    | The DICOM device can store (archive) DICOM data.  | <input type="radio"/>    |
| Move                     | The DICOM device can be asked to transmit DICOM data to other devices (third-party move).   | <input type="radio"/>    |
| Storage Commit           | The DICOM device is a server that supports Storage Commit, a feature that ensures that the server properly receives each transmission.  | <input type="radio"/>    |
| Query                    | When DICOM Host is [Yes], the default, and Query is checked, the DICOM device is a server that supports standard DICOM queries. When DICOM Host is [No], and Query is checked, the DICOM device is a server that supports SQL queries.                      | <input type="radio"/>    |
| Relational               | The DICOM device is a server that supports relational query types. When not selected, only hierarchical queries are supported.  | <input type="radio"/>    |
| PPS Manager              | (Available only when IHE Support is enabled. This option must be selected for exactly one host.)<br>Select this option to designate this host as a Performed Procedure Step Manager. This activates redistribution of Performed Procedure Step information. | <input type="radio"/>    |
| Security TLS             | Not used. Do not check the box.   | <input type="checkbox"/> |

| Item              | Description  | Mod |
|-------------------|--|-----|
| Transfer Syntaxes | <p>(Available only for hosts with Store capability.)</p> <p>Configured separately on the Acquisition and Reference computers, defines which transfer syntaxes are to be negotiated for Send operations from the Acquisition and Reference computers to this host. The syntaxes do not apply to Default Copy operations or the automatic Acquisition to Reference transfers.</p> <p>Click [Transfer Syntaxes]. The Transfer Syntaxes Supported dialog box appears.</p>  <p>There are three syntax groups, [XA] for loops, [SECONDARY] for still images (Reference and Annotation), and SR for RDSR. Within each group, there are two uncompressed syntaxes, [LEI] for little endian implicit, and [LEE] for little endian explicit, plus a [JPEG] syntax.</p> <p>Select all desired syntaxes in both the [XA] and [SECONDARY] groups (check marks appear).</p> <p>Select all desired syntaxes in both the [XA] and [SECONDARY] groups (check marks appear).</p> <p>Then, with the mouse, drag all selected syntaxes into the desired priority order (highest priority at top).</p> <p>Before the Send operation begins, all selected XA and Secondary syntaxes (that correspond to what was selected in the Studies Management window) are proposed to the host. Then, within each group (XA and Secondary), the highest-priority syntax that is accepted by the host, [that is already in the cache], is used for the transfer.</p> <p>Even if it is not selected, an implied LEI selection (lowest priority) always exists.</p> <p>If one of the selected syntaxes is also the native syntax (typically Lossless) then it is given priority over all other selected syntaxes.</p> <p>If the host rejects all proposed syntaxes, then LEI is used. If in the unusual circumstance, LEI is rejected, an error occurs and nothing is sent.</p> <p>If no syntaxes are selected, then whatever was selected in the Studies Management window (first priority) plus LEI are proposed to the host. If the Studies Management selection is rejected, LEI is used.</p> | ○   |




| Item       | Description   | Mod                    |
|------------|---|------------------------|
| Processing | <div><p>Set image processing when transfer the image. Click [Processing] to display [Transfer Image Processing Configuration] window is displayed. Select an appropriate image processing and click [OK].</p><div><div><div>Transfer Image Processing Configuration</div><div>CVS_GW_01</div><div><div>Automatic Transfer</div><div><div><div><div><input checked="" type="checkbox"/> AWL</div><div><input checked="" type="checkbox"/> Edge Enhancement</div><div><input checked="" type="checkbox"/> Negative/Positive Inversion</div><div><input checked="" type="checkbox"/> Subtraction</div><div><input checked="" type="checkbox"/> Zoom</div><div><input checked="" type="checkbox"/> Resize</div><div><input type="checkbox"/> Resize DSA</div></div><div><input type="checkbox"/> Check/Uncheck All</div></div></div><div><div>Manual Transfer</div><div><div><div><input checked="" type="checkbox"/> AWL</div><div><input checked="" type="checkbox"/> Edge Enhancement</div><div><input checked="" type="checkbox"/> Flip Horizontal</div><div><input checked="" type="checkbox"/> Flip Vertical</div><div><input checked="" type="checkbox"/> Mask Weight</div><div><input checked="" type="checkbox"/> Negative/Positive Inversion</div><div><input checked="" type="checkbox"/> Pixel Shift</div><div><input checked="" type="checkbox"/> Subtraction</div><div><input checked="" type="checkbox"/> Zoom</div></div><div><input type="checkbox"/> Check/Uncheck All</div></div></div><div><div>OK</div><div>Cancel</div></div></div></div></div><div><div><div><div><div><div></div><div></div></div><div>NOTE</div></div></div><div><p>If changing the setting of auto transfer after specifying default server, the change does not reflect to ACQ PC. Specify the default server again or set the same setting to ACQ PC.</p></div></div></div></div> | <div><div></div></div> |

## DICOM Tab

DICOM options are configured as follows.

| Item                          | Description  | Mod |
|-------------------------------|--|-----|
| <b>DICOM Options</b>          |  |     |
| Hosts Authentication          | Causes this computer to use DICOM host authentication when communicating with any DICOM device.  | ×   |
| SCU Timeout                   | Defines the maximum number of seconds to wait for an answer while acting as an SCU. If communication is not successful within this maximum, an error will occur.   | ×   |
| SCP Timeout                   | Defines the maximum number of seconds to wait for an answer while acting as an SCP. If communication is not successful within this maximum, an error will occur.   | ×   |
| MWL Auto refresh delay        | Defines the number of minutes to delay between each automatic refresh of the procedure steps list when the Worklist Management dialog box is open and its Auto Refresh is enabled.   | ×   |
| MWL Skip IHE Extension Fields | Defines the description of the “scheduled procedure steps (0040,0007)” to remove from a query when acquiring patient information on the Worklist Management.   | ×   |
| <b>Storage Commit</b>         |  |     |
| Max. Requests                 | The maximum number of storage commit requests without response that can be reached for a DICOM file before generating an error. A storage commit request is considered to be without response if the response has not occurred within Timeout seconds. | ×   |
| Max. Failed                   | The maximum number of failed storage commit requests before aborting file transfer attempts and displaying an error.   | ×   |
| Interval                      | The number of seconds between each scan of studies in the local cache to verify whether there are storage commit requests to be sent. The first scan is made this number of seconds after the start of the GUI.  | ×   |
| Timeout                       | The maximum number of seconds permitted to elapse between a storage commit request being sent and a retry being made if there was no response to the request.  | ×   |

| Item   | Description  | Mod |
|--|--|-----|
| <b>RDSR (For acquisition use only)</b>       |  |     |
| Calibration Uncertainty (%)                  | Acceptable value of error from measurement value by dosimeter.   | ○   |
| Calibration Factor                           | Coefficient of dosimeter used in calibration.  | ○   |
| Calibration Responsible Party                | Name of manager or group for calibration.  | ○   |
| Calibration Protocol                         | Protocol name for calibration.   | ○   |
| Calibration Date                             | Calibration date.  | ○   |
| <b>MPPS Options</b>                          |  |     |
| Enable MPPS Support                          | Enables MPPS Support when checking the check box, and the following option can be modified.<br>If this option is selected, at least one host should be able to use PPS manager.<br> <a href="#">"DICOM Hosts Tab" P.17-20</a> | ×   |
| Server supports Unscheduled Cases            | Support MPPS for unscheduled studies when checking the check box.  | ×   |
| Include Radiation Dose module                | Check if sending dose information to MPPS server.  | ×   |
| Allow to overwrite MWL Study Description     | Activate to be able to correct [Study Description] when starting the study from MWM server with PPS manager.   | ×   |
| Use DUP Group Name for Protocol Name         | Activate to start new study without inputting protocol name on "New Study" window, when starting the study from MWM server with PPS manager.<br>In this case, default DUP group name will be used for protocol name.<br>Ex.) Cardio, Head, Abd...  | ×   |
| Nb of commands retries:                      | Maximum number of retries of IHE command if failed.  | ×   |
| Delay between retries (sec):                 | Delay between retries.   | ×   |
| Delay after fatal failure (min):             | (Applies when IHE command is not succeed though it reaches the maximum number of retries.) Defines the waiting time before repeating retries. Retries will not cancel unless the user cancels the command.   | ×   |
| Scheduled Protocol Codes:                    | Click <b>View</b> to display or edit current protocol code. Click <b>Load</b> to import the setting file and replace all protocol codes to new protocol codes of files provided from HIS.  | ×   |
| Error codes:                                 | Click <b>View</b> to display or edit current error code. Click <b>Load</b> to import the setting file and replace all error codes to new error codes of files provided from HIS. Also includes default error.  | ×   |
| <b>Presentation State Options</b>            |  |     |
| Enable Presentation State Support            | IHE support is enabled if you checked, and the following options can be used.<br>If this option is selected, at least one host must enable PPS manager.<br>See <a href="#">"DICOM Hosts Tab" P.17-20</a>   | ×   |
| Default Presentation State group:            |  | ×   |
| Create acquisition Presentation Stage group: |  | ×   |

## 17 Administration

---

| Item | Description  | Mod |
|------|--|-----|
|      | Automatically apply default Presentation State group:                                | ×   |
|      | Prompt for new Presentation State group creation:                                    | ×   |
|      | Presentation State User Mode:  | ×   |
|      | Prompt for Presentation State selection when sending a local study on the network:   | ×   |
|      | Apply acquisition Presentation State group when sending acquisition study to server: | ×   |

## Software Tab

Software options are configured as follows.

### Reference Monitor

The screenshot shows the 'Reference Monitor' configuration window. On the left is a sidebar with a tree view containing the following items: System, DICOM Hosts, DICOM, **Software** (highlighted), Hardware, Hardware (Shutter), Storage (Local), Storage (Network), Display, Menus, Notifications, Physicians, Study Information, Fonts & Colors, Database, Devices, External Software, Fusion, and Logs. The main area is titled 'Software Configuration' and contains several sections:
 

- Close Study**: Includes a 'Set default to' label and two checked checkboxes: 'Close study in Acquisition' and 'Close study in Review'.
- Print Options**: Includes two checked checkboxes: 'Print in Color (Non-DICOM)' and 'Keep the Print Manager open after printing'.
- Image Selector Filtering (from C-Arm position)**: Includes a 'Degree of tolerance' input field with the value '5'.
- View Direction**: Includes a 'Frontal' label and two dropdown menus, both set to 'Detector To Source'.
- Automatic Plane Mode Change**: Includes a checked checkbox labeled 'Enable'.
- Remote Maintenance**: Includes a checked checkbox 'Site-View Plus', a 'URL' input field with 'http://192.168.100.100/jp', a checked checkbox 'MPC Error Report', and an 'MPC' input field with '192.168.100.100'.




 At the bottom right are three buttons: 'OK', 'Cancel', and 'Apply'.


### Acquisition Monitor

The screenshot shows the 'Acquisition Monitor' configuration window. It has the same sidebar as the Reference Monitor window. The main area is titled 'Software Configuration' and contains:
 

- Image Selector Filtering (from C-Arm position)**: Includes a 'Degree of tolerance' input field with the value '5'.
- View Direction**: Includes a dropdown menu set to 'Detector To Source'.
- Automatic Plane Mode Change**: Includes an unchecked checkbox labeled 'Enable'.
- Sounds**: A button located to the right of the View Direction section.

 At the bottom right are three buttons: 'OK', 'Cancel', and 'Apply'.

| Item                               | Description  | Mod |
|------------------------------------|--|-----|
| <b>Close Study</b>                 |  |     |
| Close Study                        | (Reference only) Defines the default settings in the Close confirmation message box that appears when a study in acquisition or review is closed. Check the options that you want to have checked by default in the confirmation box.  | ○   |
| <b>Print Options</b>               |  |     |
| Print Options                      | <p>(Reference only) Check [Print in Color] so that color elements will be printed on non-DICOM printers that support color such as laser and ink jet.</p> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;"> <div style="text-align: center;">  <b>NOTE</b> </div> <p>If you print an image that has already been saved in DICOM format, colors will appear as shades of gray.</p> </div> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;"> <div style="text-align: center;">  <b>NOTE</b> </div> <p>If you print an image that has already been saved in DICOM format, colors will appear as shades of gray.</p> </div> | ×   |
| <b>Image Selector Filtering</b>    |  |     |
| Image Select Filtering             | For when filtering the Image Selector icons by C-arm position. Degree of tolerance defines within how many degrees of the current C-arm position an image's C-arm position must be, to be included in the filtered list. Both angles must be within the tolerance. For example, if Degree of tolerance is <b>5</b> and the C-arm is at position LAU-35/CAU-48, images acquired at a position of LAU-39/CAU-50 will be included in the list but images from position LAU-45/CAU-55 will not.  | ○   |
| <b>View Direction</b>              |  |     |
| Frontal                            | Set eye direction from Frontal.  | ○   |
| <b>Automatic Plane Mode Change</b> |  |     |
| Automatic Plane Mode Change        | See  "10.5 Administration" P.10-15 for details.   | ×   |
| <b>Remote Maintenance</b>          |  |     |
| SiteViewPlus+                      | Not used.  | ×   |



| Item             | Description  | Mod |
|------------------|--|-----|
| MPC Error Report | <p>Set IP address of MPC.</p> <div> <b>NOTE</b></div> <p>Making a contract with remote maintenance is necessary to enable MPC Error Report.</p> | ×   |

## Hardware Tab

Hardware options are configured as follows.

| Item                 | Description  | Mod |
|----------------------|--|-----|
| <b>Information</b>   |  |     |
| Grid                 | Defines whether a grid is used in the X-ray imaging system. Select <b>IN</b> if a grid is used. This option is used to complete the DICOM information for acquired images. This option must match the installation.  | ×   |
| FPD Size             | Specifies the size (inches) of the FPD (Flat-Panel Detector) (maximum 25 inches). This option is used to complete the DICOM information for acquired images. This value must match the installation.   | ×   |
| FOD                  | Specifies the focal point to object (ISO Center) distance. Also known as SOD (Source to Object Distance). This is a fixed value (typically 72 cm) according to C-arm model. This value (cm) must match the installation. (0 to 100.)   | ×   |
| ISO Center to Skin   | The offset distance from the FOD (ISO Center) toward the focal point, to the point approximately at the patient's skin surface (top of table) at which dosage is measured. The FDA requires this value to be 15 cm.  | ×   |
| <b>Field of View</b> |  |     |
| Field of View        | Defines eight field-of-view values (inches). Specify the length of one side of each field-of-view. This option is used to complete the DICOM information for acquired images. This option must match the installation. Large and FOV 5th to FOV 8th can be no larger than FPD Size. Medium, Small, and Smallest must be progressively smaller. | ×   |



| Item                                   | Description   | Mod |
|--|---|-----|
| <b>Curves Acquisition</b>              |   |     |
| Curves Acquisition                     | <p>Defines which (check mark shown) probe signals are to be acquired at the same times as images. Channels 1 and 2 are dedicated to ECG and Blood Pressure. Channels 3 and 4 are individually selectable. Labels entered in the <b>Label</b> column appear on screen at the left edge of the probe signal (Up to 64 keyboard characters, excluding the semi-colon ";").</p> <p>[Default Acquire Waves] makes [Acquire Waves] selected by default for new studies. See  <a href="#">"4.3.2 Defining a New Study" P.4-10.</a></p> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;"> <div style="text-align: center;">  <b>NOTE</b> </div> <p>For meaningful information, make sure that what is configured here matches the hardware installation.</p> </div> | ×   |
| <b>Dose Meter Configuration</b>        |   |     |
| Dose Meter Configuration               | Select a dosimeter connected to DAR-9500f.  | ×   |
| <b>Table Configuration</b>             |   |     |
| Table Reposition                       | If check this box, table reposition becomes available. Check if Catheterization Table KS-100 is installed.  | ×   |
| Show ROI indicating the table position | When put a check mark, Positioning mode becomes active as a default.  | ×   |
| Stroke Length (mm)                     | Set the length of stroke in the longitudinal direction of the table.  | ×   |
| <b>Large Monitor Configuration</b>     |   |     |
| Large Monitor Control                  | Check this box if SMART Display is installed.   | ×   |
| LMM IP                                 | Set IP address of SMART Display.  | ×   |
| <b>SMART Touch Configuration</b>       |   |     |
| SMART Touch Configuration              | Set unit name of SMART Touch. Select a connected unit.  | ×   |

## Hardware (Shutter) Tab


Hardware (Shutter) options are configured as follows.



|          | Frontal   |       |           |       |           |       |
|----------|-----------|-------|-----------|-------|-----------|-------|
|          | C-Shutter |       | H-Shutter |       | V-Shutter |       |
|          | Open      | Close | Open      | Close | Open      | Close |
| Large    | -15       | 62    | -15       | 100   | -15       | 100   |
| Medium   | -16       | 56    | -15       | 100   | -15       | 100   |
| Small    | -18       | 50    | -15       | 100   | -15       | 100   |
| Smallest | -24       | 35    | -15       | 100   | -15       | 100   |
| FOV 5th  | 0         | 100   | 0         | 100   | 0         | 100   |
| FOV 6th  | 0         | 100   | 0         | 100   | 0         | 100   |
| FOV 7th  | 0         | 100   | 0         | 100   | 0         | 100   |
| FOV 8th  | 0         | 100   | 0         | 100   | 0         | 100   |


| Item              | Description  | Mod |
|-------------------|--|-----|
| Shutter Positions |  |     |
| Shutter Positions | Defines the cross-direction collimator (C-shutters), horizontal (H-shutters) and vertical (V-shutters) shutters position range for each field of view of the installation. Position values are normalized to 0 to 100, where smaller values represent more open shutters (more x-ray exposure) and larger values represent more closed shutters (less X-ray exposure). [Open] must be set to a value lower than [Close]. | ×   |

## Storage (Local) Tab

Storage (Local) options are configured as follows.




| Item                     | Description  | Mod |
|--------------------------|--|-----|
| <b>Cache Information</b> | Hard disk space is used for both data cache (folder \Cache) and network cache (\Cache\NetCache). Incorrect cache settings can cause unexpected system behavior so DO NOT modify items in this box.   |     |
| Open                     | Open cache folder.   | ×   |
| Cache Selector           | <p>Choose the letter of the hard drive on which the two cache folders will reside. Folders \Cache and \Cache\NetCache are automatically created on the chosen drive.</p> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;">  <b>NOTE</b><br/>           The cache must reside on a separate dedicated drive or drive array. Other software must never access this drive.         </div> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;">  <b>NOTE</b><br/>           If you change the drive letter for the local cache and wish to preserve existing studies, exit the software, move all data from the old \Cache folder to the new one, and then restart the software.         </div> | ×   |

| Item                         | Description  | Mod |
|------------------------------|--|-----|
| Lower and Upper Delete Limit | <p>[Upper delete limit] defines the minimum amount of space you want to keep free on your Cache drive. Whenever free space drops below this threshold, unprotected studies are deleted, starting with the least-recently accessed, until there is at least this amount of free space on the Cache drive.</p> <p>When a study is open for acquisition, deletion of studies only starts when the amount of free space falls below [Lower delete limit].</p> <p>Both these values are expressed in MB, where 1MB equals 1024x1024 or 1,048,576 bytes. For each GB of cache space, enter 1024 MB. For example, for 50 GB, enter 51200 MB (50*1024).</p> <p>The difference between [Upper delete limit] and [Lower delete limit] should be a little larger than the largest study that you ever expect to create. Here are guidelines for setting the [Upper delete limit] assuming the following:</p> <ul style="list-style-type: none"> <li>• MS=Maximum size of one study (e.g., 2000 MB)</li> <li>• MN=Maximum space consumed in NetCache (e.g., 3000 MB)</li> <li>• DF=Desired free space remaining with full cache (e.g., 5000 MB)</li> </ul> <p>[Upper delete limit] is calculated as MS+MN+DF, or 10000 MB.</p> | ×   |
| <b>Primary Archive</b>       |  |     |
| Manual Archive               | Not used.  | ×   |
| Primary Archive              | Designates the local system as the primary archive system for your department or institution. In this mode, the system must write a study to CD/DVD before it can be unprotected and subsequently deleted.   | ×   |
| CD/DVD Copies                | (Primary Archive mode only) Defines how many times the local system must write a study to CD/DVD before it is automatically unprotected and subject to deletion. (1 to 9.)   | ○   |
| Cache Warning Threshold (MB) | <p>(Primary Archive mode only)</p> <p>This value, specifies in megabytes, (default 8000) the minimum amount of free space required on the local system cache drive to avoid having a low-cache warning being entered in the DAR-9500f message log (question-mark icon in status bar), and if configured, a warning message being sent to an administrator's computer. Do not enter a value of less than the [Upper delete limit] value of the Cache Information (System tab). The Cache Full Shell Command can be used to display a warning message on the administrator's monitor.</p> <div style="border: 1px solid black; padding: 10px; margin-top: 10px;"> <div style="text-align: center;">  <b>NOTE</b> </div> <p>To free up space, you can delete studies from the local system (Studies Management window) that have already been written to CD/DVD the required number of times. You can also write studies to CD/DVD and then delete them from the local system.</p> </div>  | ○   |
| <b>File Management</b>       |  |     |
| File Management              | (Reference only) Enables [Allow individual file transfer] function, when checking the check box.   | ○   |


## Storage (Network) Tab


Storage (Network) options are configured as follows.

Reference Monitor



Acquisition Monitor




| Item  | Description  | Mod                      |
|---|--|--------------------------|
| <b>Default Server</b>   |  |                          |
| Images  | (When <b>Primary Archive</b> mode not selected) Select the server from the Default Server drop-down list. All newly-acquired images during studies and the dose report images created after study will be sent to the archive server. Choose this <b>ONLY</b> on the Reference computer.   | <input type="radio"/>    |
| RDSR  | (When <b>Primary Archive</b> mode not selected) Newly-created RDSR reports will be sent to the selected server. Choose this <b>ONLY</b> on the Reference computer.   | <input type="radio"/>    |
| <b>Send Format</b>  |  |                          |
| Send Format   | (Acquisition only) Usually check lossless.   | <input type="checkbox"/> |
| <b>Default Copy Destination</b>                               |  |                          |
| Default Copy Destination                                      | <p>(Available only if the default server supports "Move" operations.) The check box list shows every station with Store capability <b>except the default server</b> as defined in DICOM Hosts AE Table of the Hosts tab. Check each station in the list to which you want a copy of the study sent as soon as the study in acquisition is closed. This is in addition to the study being first transferred from the Acquisition station to the default server. The completed study is copied from the default server to every station checked in this list.</p> <p>Furthermore, you can define which study series types are to be included in the copy. To do this, right-click a station in the list. The Series pop-up appears. Check only the series types you want sent and then click <b>Close</b>. For example, you could choose to send only annotated and reference images to a reporting server.</p>  | <input type="checkbox"/> |
| <b>3D Workstation</b>   |  |                          |
| 3D Workstation  | Choose the name of the 3D workstation to which each new acquisition will be sent.  | <input type="checkbox"/> |
| <b>Processed Image Transfer</b>                               |  |                          |
| Check "Process Before Transfer" in Manual Transfer as default | Transfer the image after processing when sending individually.   |                          |
| Send processed images only                                    | Only processed images will be transferred.   |                          |
| Exclude DSA from Process Before Auto Transfer                 | Transfer the DSA image without image processing for auto transfer.   |                          |
| Force Process Before Auto Transfer of StentView/StentShot     | Transfer the image by forcing the image processing on StentView and StentShot image for auto transfer.   |                          |


■ Display Tab


Display options are configured as follows.

Reference Monitor



Acquisition Monitor



| Item                                   | Description   | Mod                                 |
|--|---|-------------------------------------|
| <b>Display Properties</b>              |   |                                     |
| Default to Play All Loops              | (Reference only) Causes loops to be played back one after the other instead of just playing the selected loop continuously.   | <input type="radio"/>               |
| Display Image-Selector Icons           | Causes the Image-Selector icons to be displayed on the monitor right edge, either in a single column, the default, or in a double column when [Double Rows or Columns] is checked. Leave [Vertical] checked. C-arm position is displayed on the image selector when [Show C-Arm Position] is checked.   | <input checked="" type="checkbox"/> |
| Display Information                    | Displays the study information line at the top left of the Image Viewer window, optionally in a Large Font.   | <input checked="" type="checkbox"/> |
| Synchronize on Loops                   | (Reference only) If the Reference computer is displaying the study currently in acquisition (loops or stills), when a new acquisition is received by the Reference computer, it will begin playing the just-acquired loop.  | <input type="radio"/>               |
| Default REF Sync Mode                  | (Reference only) Frontal Reference monitor and Lateral Reference monitor are synchronized.  | <input type="radio"/>               |
| Default ACQ Sync Mode                  | ACQ Frontal and ACQ Lateral are synchronized.   | <input type="radio"/>               |
| Automatically Display Analysis Results | (Reference only) Displays the QCA and LV analysis results in the upper left of the Image Viewer window.   | <input type="radio"/>               |
| Display Active Status                  | (Acquisition only) (Off by default) Displays the floating status icon in the Image Viewer window. This has been replaced by the Acquisition Active message that is displayed in the upper-right area of the Image Viewer window during Fluoro and Rad acquisition.  | <input type="radio"/>               |
| Monochrome Monitor                     | Optimize the side menu GUI buttons for monochrome monitor. Normally check this box.   | <input checked="" type="checkbox"/> |
| PCI View mode                          | If check this box, display loops selected before X-ray exposure.  | <input type="radio"/>               |
| Show Low Accuracy C-Arm Position       | Acquisition with Precession and Pendulum C-arm motion can also save angle information for each frame. However, there is an error up to 3 degrees between them.  | <input type="radio"/>               |
| <b>Virtual Collimator</b>              |   |                                     |
| Virtual Collimator                     | <p>Virtual collimators are displayed on the Acquisition monitor showing the position of the hardware collimators.<br/>           [Delay (sec)] defines the time in seconds (0 to 100) that virtual collimators are displayed after the hardware collimators have stopped moving.<br/>           [Width] defines the thickness of the displayed line (range 2 to 30).</p> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;"> <div style="display: flex; align-items: center; justify-content: center;">  <div style="margin-left: 10px;"><b>NOTE</b></div> </div> <p>Virtual collimation function indicates the collimator leaf position before fluoroscopy as a suggestion. Therefore, the position may misalign by the field of view size, C-arm position, and etc.</p> </div> | <input type="radio"/>               |




| Item              | Description   | Mod |
|-------------------|---|-----|
| <b>Display In</b> |   |     |
| Display In        | <p>Choose a parameter in the list and then check any available combination of Fluoro, Rad, Playback, Review or Displayed. Some parameters can only be displayed at certain times, as follows:</p> <ul style="list-style-type: none"> <li>• Rad: (Acquisition monitor only) During actual Rad acquisition. Applicable only when a study is open for acquisition.</li> <li>• Fluoro: (Acquisition monitor only) During actual Fluoro acquisition, and if configure, during hold of last Fluoro image. Applicable at all times.</li> <li>• Playback: (Acquisition monitor only) When a Rad loop or an optionally-recorded Fluoro loop is automatically played back, or the last Fluoro image frame is being displayed, both after an acquisition. Applicable only when a study is open for acquisition.</li> <li>• Review: For all playback on Reference monitor. (Acquisition monitor only) When both a study is open for acquisition and a loop is manually selected for playback.</li> <li>• Displayed: For Study dose, AK (Acquisition monitor only) and Storage (Acquisition monitor only).</li> </ul> <p>The correspondence between parameter names in this list and their on-screen labels are as follows:</p> <p style="padding-left: 40px;">dGym<sup>2</sup>: Study Dose (Total dosage since start of study),<br/>(Acquisition monitor) and Loop Dose (Dose of the loop), (Reference monitor)</p> <p style="padding-left: 40px;">mGy/min: AKR</p> <p style="padding-left: 40px;">kV: kV</p> <p style="padding-left: 40px;">mA: mA</p> <p style="padding-left: 40px;">ms: Exp Time</p> <p style="padding-left: 40px;">LAO or RAO: LAO-RAO</p> <p style="padding-left: 40px;">CRA or CAU: CRA-CAU</p> <p style="padding-left: 40px;">SID (cm): SID</p> <p style="padding-left: 40px;">IA Delay (s): IA-Delay</p> <p style="padding-left: 40px;">FOV (in): FOV</p> <p style="padding-left: 40px;">Fluoro (min): Fluoro Time</p> <p style="padding-left: 40px;">Acq Inj Time (s): Acq Inj Time</p> <p style="padding-left: 40px;">Storage (min): Storage</p> <p style="padding-left: 40px;">&lt;custom&gt;: Channel 1 (probe signal)</p> <p style="padding-left: 40px;">&lt;custom&gt;: Channel 2 (probe signal)</p> <p style="padding-left: 40px;">&lt;custom&gt;: Channel 3 (probe signal)</p> <p style="padding-left: 40px;">&lt;custom&gt;: Channel 4 (probe signal)</p> <p style="padding-left: 40px;">&lt;text&gt;: Image Comment</p> <p style="padding-left: 40px;">&lt;seconds&gt;: Fluoro Map Time</p> <p style="padding-left: 40px;">Fram Inj Time(s): Frame Inj Time</p> <p style="padding-left: 40px;">Acquisition Time: Acq Time</p> <p style="padding-left: 40px;">&lt;LIVE/LIH/ STORED&gt;: Image Displays</p> <p style="padding-left: 40px;">&lt;L/R/A/P&gt;: Orientation 1</p> <p style="padding-left: 40px;">&lt;F/H/A/P&gt;: Orientation 2</p> | ○   |


| Item                           | Description  | Mod                   |
|--------------------------------|--|-----------------------|
| Study Information Display      |  |                       |
| Accession Number               | Displays at the very top of Study Information on the image viewer at the start of new study.                   | <input type="radio"/> |
| Study ID                       | Displays at the very top of Study Information on the image viewer at the start of new study.                   | <input type="radio"/> |
| Patient ID                     | Displays at the very top of Study Information on the image viewer at the start of new study.                   | <input type="radio"/> |
| Patient Name                   | Displays at the very top of Study Information on the image viewer at the start of new study.                   | <input type="radio"/> |
| Study Date                     | Displays at the very top of Study Information on the image viewer at the start of new study.                   | <input type="radio"/> |
| Patient Age                    | Displays at the very top of Study Information on the image viewer at the start of new study.                   | <input type="radio"/> |
| Patient Sex                    | Displays at the very top of Study Information on the image viewer at the start of new study.                   | <input type="radio"/> |
| Patient Birth Date             | Displays at the very top of Study Information on the image viewer at the start of new study.                   | <input type="radio"/> |
| Image Information to Display   |  |                       |
| Image Type                     | Displays at the upper-right of Study Information on the image viewer.  | <input type="radio"/> |
| Total Frames                   | Displays at the upper-right of Study Information on the image viewer.  | <input type="radio"/> |
| Total Loops                    | Displays at the upper-right of Study Information on the image viewer.  | <input type="radio"/> |
| IVR NEO Status                 |  |                       |
| Review (Reference only)        | Causes the indicator “Ref” to appear in upper-left corner of each monitor that is in Reference mode (IVR NEO). | <input type="radio"/> |
| Acquisition (Acquisition only) |  | <input type="radio"/> |
| DSA Tools                      |  |                       |
| Mask Weight Increment          | Set an increment for the Up/Down button of "Mask Weight" used in DSA tools. (1 to 25)<br>Default value is "1". | <input type="radio"/> |

## Menus Tab

Menu options are configured as follows.




| Item         | Description   | Mod                   |
|--------------|---|-----------------------|
| Side Menu    | Causes the side menu at the left edge of the Image Viewer window to be automatically hidden after [Delay (sec.)] seconds (1 to 10) without mouse activity. Moving the mouse to the left screen edge restores both the side menu and bottom status bar.      | <input type="radio"/> |
| Bottom Menu  | Causes the Cine Control / Status bar along the bottom of the Image Viewer window to be automatically hidden after [Delay (sec.)] seconds (1 to 10) without mouse activity. Moving the mouse to the bottom screen edge, restores only the bottom status bar. | <input type="radio"/> |
| Menu Visible | Causes a thin gray bar to be shown in place of the hidden side menu and / or bottom status bar, indicating where you can move the pointer to restore the hidden elements.   | <input type="radio"/> |

For information about items in the [Menu Creation] group, see  ["17.7.4 Menus and DUP Configuration" P.17-69](#).

## Notifications Tab

Notifications options are configured separately on the Reference and Acquisition computers as follows.

### Reference Monitor




The Reference Monitor Notifications Configuration window is shown. The left sidebar lists various system settings, with 'Notifications' selected. The main area is titled 'Notifications Configuration' and contains three sections:

- Display Message / Warning for:**
  - ☐ Print list image addition
  - ☒ Reset calibration factor
  - ☒ Show the calibration range
  - ☒ Study importation completed
- Display Confirmation for:**
  - ☒ Analysis calibration factor retention
  - ☒ Delete created objects
  - ☒ Print list all image deletion
  - ☒ Print list single image deletion
  - ☒ Protected study overwrite
  - ☒ Retrieve study data from server before writing media
  - ☒ Study importation
  - ☒ Study Management list refresh after study deletion
  - ☒ Study visualization
- Pediatrics:**
  - ☐ Show pediatric warning
  - Upper age limit: 21 years
- Dose:**
  - ☒ Show dose warning
  - Warning Threshold: 2000 mGy

At the bottom right are buttons for 'OK', 'Cancel', and 'Apply'.

### Acquisition Monitor



The Acquisition Monitor Notifications Configuration window is shown. The left sidebar lists various system settings, with 'Notifications' selected. The main area is titled 'Notifications Configuration' and contains three sections:

- Display Message / Warning for:**
  - ☒ Do not move catheter
  - ☒ Not ready to record (enter patient information)
- Display Confirmation for:**
  - ☒ Protected study overwrite
  - ☒ Study Management list refresh after study deletion
  - ☒ Study visualization
- Pediatrics:**
  - ☐ Show pediatric warning
  - Upper age limit: 21 years
- Dose:**
  - ☐ Show dose warning
  - Warning Threshold: 2000 mGy

At the bottom right are buttons for 'OK', 'Cancel', and 'Apply'.

Select the messages that you wish to have shown in prompt boxes. If you do not select a message, the user can usually rely on equivalent messages displayed in the status bar.

| Item                   | Description  | Mod                   |
|------------------------|--|-----------------------|
| <b>Pixel Spacing</b>   |  |                       |
| Pixel Spacing          | Set acceptable value of Pixel Spacing. Setting percentage is the acceptable value for calibration result.        | <input type="radio"/> |
| <b>Pediatrics</b>      |  |                       |
| Show pediatric warning | If check this box, warning is displayed at the start of new study if patient age is under the [Upper age limit]. | <input type="radio"/> |
| <b>Dose</b>            |  |                       |
| Show dose warning      | If check this box, warning is displayed on the Acquisition monitor if dose exceeds the [Warning Threshold].      | <input type="radio"/> |

## Physicians Tab



[Physicians] configuration is performed only on the Reference computer.

The [Physicians] tab enables you configure options specific to each performing physician, including the default. For each physician, you can select the default Radiography and Fluoroscopy programs and Fluoro record mode, and configure the IVR NEO/SMART Touch. Before setting options here, first enter all performing physicians by adding them to the [Performing Physicians Table] found on the [Database] tab. Physician options are configured as follows.

Selecting Shortcuts keys Configuration:

The screenshot shows the 'Configuration' window with the 'Physicians Configuration' tab selected. The 'Shortcuts keys' section is active, showing a list of keys (F1-F8) and their corresponding functions. The 'F1' key is selected, and its function is 'Select FLUORO Program'. The 'Parameters' section is empty. The 'Assign' and 'Delete' buttons are visible. The 'Static Threshold (%)' is set to 15, and the 'Step/Play Threshold (%)' is set to 30. The 'IVR NEO' and 'SMART Touch' settings are also visible.

Selecting SMART Touch Configuration:

The screenshot shows the 'Configuration' window with the 'Physicians Configuration' tab selected. The 'SMART Touch' section is active, showing a list of keys (Panel B1-B8) and their corresponding functions. The 'Panel B1' key is selected, and its function is 'Select FLUORO Program'. The 'Parameters' section is empty. The 'Assign' and 'Delete' buttons are visible. The 'Unit' is set to 'TOUCH-01-01', and the 'Layout' is set to '5X4'. The 'Color' is set to 'Dark'. The 'Copy current touch config for Default / TOUCH-01-01 to' and 'Copy Destination Physician' and 'Copy Destination Unit' fields are also visible.

Selecting Road MAP Fluoro Configuration:

**Configuration**

**Physicians Configuration**

Physician Name:

E-Shutters: ☐ Display E-Shutters, Opacity (%):

Study Loops: ☒ Auto Transfer

Acquisition Default Menu: ☒ Radiography, ☐ Fluoroscopy

Fluoro Record: ☐ Direct, ☒ Last N Seconds

Acquisition: ☒ Hold Fluoro, ☐ Hold Recorded Fluoro

Shortcut keys: ☒ Coordination with Geometry

Enhance Contour: Region Rendering Quality:  Smoothing:  Thickness: ☐ Thin, ☒ Medium, ☐ Thick

Mask Region: Initial Size (%):  Minimum Size (%):

Sketch Edition: Dots Size:  Line Type:  Line Width:  Line Color:  Default Line:

OK Cancel Apply

Selecting Geometry Configuration:

**Configuration**

**Physicians Configuration**

Physician Name:

E-Shutters: ☐ Display E-Shutters, Opacity (%):

Study Loops: ☒ Auto Transfer

Acquisition Default Menu: ☒ Radiography, ☐ Fluoroscopy

Fluoro Record: ☐ Direct, ☒ Last N Seconds



Acquisition: ☒ Hold Fluoro, ☐ Hold Recorded Fluoro

Shortcut keys: ☒ Coordination with Geometry




Sketch Display Timeout (sec):  Oblique Angle Margin (degree):  Sagittal Angle Margin (degree):  SID Margin (mm):


Table Lateral Position Margin (mm):  Table Longitudinal Position Margin (mm):  Table Height Margin (mm):  Ceiling Travel Lateral Position Margin (mm):  Ceiling Travel Longitudinal Position Margin (mm):  Table Head Tilt Angle Margin (degree):  Table Cradle Tilt Angle Margin (degree):  Table Horizontal Rotation Angle (degree):

OK Cancel Apply

| Item                     | Description   | Mod                   |
|--------------------------|---|-----------------------|
| Physician Name           | Choose a physician from [Physician Name].   | <input type="radio"/> |
| Acquisition Default Menu | <ul style="list-style-type: none"> <li>Set [Study Loops]. To set the default [Auto Transfer] setting for new studies, select [Auto Transfer]. When selected, study loops are automatically sent to the default server.</li> <li>To set the default Fluoro program, select [Fluoroscopy], click [Select], navigate through the displayed program menu and click the desired default program (selected item will be highlighted).</li> <li>To see what programs are set for a particular physician, choose the physician from the [Physician] list, and then click [Select]. The selected program is indicated by the highlighted program menu item.</li> </ul> | <input type="radio"/> |
| Fluoro Record            | Set the default Fluoro Record mode for the physician, either [Direct] or [Last N Seconds]. [Direct], requires the physician to choose when to record. [Last N Seconds] causes Fluoro images to always be recorded, allowing the physician to decide when to [Save].   | <input type="radio"/> |
| Acquisition              | Select [Hold Fluoro] to have the last Fluoro image acquired held on the Acquisition monitor after each unrecorded Fluoro acquisition. If [Hold Fluoro] is not enabled, the Acquisition monitor reverts to whatever was displayed before the Fluoro acquisition (for example, a playing Rad loop). (Only applicable to unrecorded Fluoro when a study is open for acquisition.)  | <input type="radio"/> |
| E-Shutters               | Set [E-Shutters].   | <input type="radio"/> |
| Study Loops              | Set [Study Loops]. To set the default [Auto Transfer] setting for new studies, select [Auto Transfer]. When selected, study loops are automatically sent to the default server. See  <a href="#">"4.3.2 Defining a New Study" P.4-10.</a>  | <input type="radio"/> |
| <b>Shortcuts keys</b>    |   |                       |
| Description              | Indicates the explanations of functions for selected keys.  |                       |
| Keys                     | Selects function keys and IVR NEO buttons to be assigned. A letter "F" indicates the function key and "B" for IVR NEO button. To remove the function assigned a key, select the key and click the [Delete] button.  | <input type="radio"/> |
| Functions                | Selects functions to be assigned for selected keys. See  <a href="#">"4.8.1 IVR NEO Buttons" P.4-66.</a>   | <input type="radio"/> |
| Parameters               | Enables to assign some functions with parameters if necessary. Selects parameters and click [Apply] button.   | <input type="radio"/> |
| Static Threshold         | Defines percent of maximum joystick deflection (from center position) that is permitted to occur without any action being taken. This allows for small unintentional joystick movements around the center position to be ignored. Entered as percent from 0 to 100. This value must always be smaller than [Step/Play Threshold].   | <input type="radio"/> |
| Step/Play Threshold      | Defines percent of maximum joystick deflection (from center position to left or right) to be used as the threshold for where joystick movement changes from single frame forward / backward, to variable-rate frame-by-frame mode. Entered as percent from 0 to 100. This value must always be larger than [Static Threshold].  | <input type="radio"/> |



| Item                       |                          | Description  | Mod                                 |
|----------------------------|--------------------------|--|-------------------------------------|
| Mouse Speed                |                          | Optionally adjust joystick [Mouse Speed]. This option defines how fast on the on-screen pointer is moved with the joystick when in mouse mode, for example, during angiographic analysis (QCA/LVA). Entered as a number from 2 to 100.   | <input type="radio"/>               |
| Autorepeat Time            |                          | Optionally adjust joystick [Autorepeat Time]. This option defines the amount of time (milliseconds) to wait before moving to the previous / next icon in the Image Selector when holding the joystick in the up or down position. The default of 500 ms means that when the joystick is continuously held up or down, each loop will be played for 0.5 seconds before skipping to the next one. Entered as a number from 100 to 2000 ms. | <input type="radio"/>               |
| Joystick Maximum Volume    |                          | Set Joystick Maximum Volume.<br><div style="border: 1px solid black; padding: 10px; margin: 10px 0; text-align: center;">  <b>NOTE</b><br/><br/>           Do not change the [Joystick Maximum Volume] (range 1 to 100) as this is preset according to the hardware-specific value for the joystick model installed.         </div>                     | <input type="radio"/>               |
| <b>SMART Touch</b>         |                          |  |                                     |
| Key                        |                          | Selects a function and parameter associated with the function (if any) and click [Assign] button. Repeat this process for each button. To remove the function assigned a button, select the button and click the [Delete] button.  | <input type="radio"/>               |
| Unit                       |                          | Specifies SMART Touch which is going to be applied the displayed setting.<br>DAR-9500f can connect with several SMART Touch units. Setting of each SMART Touch is accessible from REF-PC.  | <input type="radio"/>               |
| Layout                     |                          | Specifies the number of buttons on Custom Panel. The default is "5x4" and 20 buttons are displayed.  | <input type="radio"/>               |
| Color                      |                          | Specifies the color style of SMART Touch application. The default is Dark.   | <input type="radio"/>               |
| Physician Name             |                          | Specifies the physician being copied to.   | <input type="radio"/>               |
| Unit                       |                          | Specifies the SMART Touch unit as the destination of the setting. Click [Go].  | <input type="radio"/>               |
| <b>Road MAP Fluoro</b>     |                          |  |                                     |
| Coordination with Geometry |                          | When checked, SIMAP Mask image and the Sketch will hide.   | <input checked="" type="checkbox"/> |
| Contour Enhancement        | Region Rendering Quality | Specify the default value of contour range. See  <a href="#">"9.5 Contour Enhancement" P.9-8.</a>   | <input type="radio"/>               |
|                            | Smoothing                | The setting value will be used as default value. See  <a href="#">"9.5 Contour Enhancement" P.9-8.</a>  | <input type="radio"/>               |
|                            | Thickness                | Specify default value of contour thickness.  | <input type="radio"/>               |

| Item   |                  | Description   | Mod |
|--|------------------|---|-----|
| Sketch Edition                                   | Dots Size        | Specify dots size when placing the sketch.  | ○   |
|  | Line Type        | Specify line style.   | ○   |
|  | Line Width       | Specify line width.   | ○   |
|  | Line Color       | Specify line color.   | ○   |
|  | Default Line     | Specify the type which is selected when [Edit] window is displayed.   | ○   |
| Mask Region                                      | Initial Size (%) | Modify the initial size when setting ROI. Initial size is specified by percentage of the image. If the initial value is 30, the size of ROI will be displayed as 30% in vertically and horizontally. See  "9.6 Mask Region" P.9-11 | ○   |
|  | Minimum Size (%) | Set minimum size of ROI.  | ×   |
| <b>Geometry</b>                                  |                  |   |     |
| Sketch Display Timeout (sec)                     |                  | If communication is disconnected, specify the time (sec) before invalid the sketch automatically. If input "0", it will not invalid though the communication is disconnected.   | ×   |
| Oblique Angle Margin (degree)                    |                  | If Coordination with Geometry is checked and the C-arm is changed more than input angle in an oblique direction, SIMAP and Sketch Line will hide.   | ×   |
| Sagittal Angle Margin (degree)                   |                  | If Coordination with Geometry is checked and the C-arm is changed more than input angle in a sagittal direction, SIMAP and Sketch Line will hide.   | ×   |
| SID Margin (mm)                                  |                  | If Coordination with Geometry is checked and SID is changed more than input value, SIMAP and Sketch Line will hide.   | ×   |
| Table Lateral Position Margin (mm)               |                  | If Coordination with Geometry is checked and the table moves more than input value in a transversal direction, SIMAP and Sketch Line will hide.   | ×   |
| Table Longitudinal Position Margin (mm)          |                  | If Coordination with Geometry is checked and the table moves more than input value in a longitudinal direction, SIMAP and Sketch Line will hide.  | ×   |
| Table Height Margin (mm)                         |                  | If Coordination with Geometry is checked and the table height is changed more than input value, SIMAP and Sketch Line will hide.  | ×   |
| Ceiling Travel Longitudinal Position Margin (mm) |                  | If Coordination with Geometry is checked and the C-arm moves more than input value in a longitudinal direction, SIMAP and Sketch Line will hide.  | ×   |
| Ceiling Travel Transversal Position Margin (mm)  |                  | If Coordination with Geometry is checked and the C-arm moves more than input value in a transversal direction, SIMAP and Sketch Line will hide.   | ×   |
| Table Head Tilt Angle Margin (degree)            |                  | (If Catheterization Table KS-100 is installed)<br>If Coordination with Geometry is checked and the table moves more than input value in a tilt direction, hide SIMAP and Sketch Line.   | ×   |
| Table Cradle Tilt Angle Margin (degree)          |                  | (If Catheterization Table KS-100 is installed)<br>If Coordination with Geometry is checked and the table moves more than input value in a cradle direction, hide SIMAP and Sketch Line.   | ×   |
| Table Horizontal Rotation Angle (degree)         |                  | (If Catheterization Table KS-100 is installed)<br>If Coordination with Geometry is checked and the table moves more than input value in a rotation direction, hide SIMAP and Sketch Line.   | ×   |



Do not re-assign functions for the B keys. These keys are associated to the IVR NEO for which button assignments are preset.



The default physician options are used if a performing physician is not selected when creating a new study.

## Study Information Tab

Study Information options are configured as follows.



Study Information is performed only on the Reference computer.

| Item                                   | Description   | Mod                   |
|--|---|-----------------------|
| <b>Mandatory Fields</b>                |   |                       |
| Patient Information, Study Information | Defines which fields in the New Study dialog box are to be made mandatory. Set a check mark to indicate <i>mandatory</i> or clear a check mark to indicate <i>optional</i> .  | <input type="radio"/> |
| <b>Accession Number</b>                |   |                       |
| Prefix                                 | Defines the 6-character (maximum) alphanumeric value used as the prefix for the auto-generated accession number that will be used to identify each study loop. If the [Site] field is used, the value entered there is copied into this field. (Up to 6 characters: alpha, number, underscore "_", period ".")  | <input type="radio"/> |
| Current                                | (Enter only on Reference) Shows and optionally sets the current auto-incremented portion of the auto-generated Accession number. If [Prefix] is used, it appears in the accession number before this generated Current value. Accession numbers, including prefix, can be up to 16 digits in length. Make sure that you do not enter a value that will cause an Accession number to not be unique. For example, do not reset this number to a lower value that was used before. (0 to 10,000,000 (10 million).) | <input type="radio"/> |
| Incr.                                  | (Enter only on Reference) Sets the increment that is added to the [Current] value when automatically generating the next Accession Number. Alternatively, set this value to zero to turn off automatic Accession Number generation, so, for example, it can be set by an HIS/RIS system, or manually. (0 to 100.)   | <input type="radio"/> |

| Item   | Description  | Mod                      |
|--|--|--------------------------|
| <b>Local Worklist</b>                          |  |                          |
| Local Worklist Enabled                         | Local Worklist Management Window is displayed when starting New Study. The study selected from Local Worklist can be deleted as the following instruction.   | <input type="radio"/>    |
| Auto Delete Study                              | Selected study will delete automatically.  | <input type="radio"/>    |
| Never Delete                                   | Selected study will not delete.  | <input type="radio"/>    |
| Prompt For Delete                              | After creating a new study, user is prompted.  | <input type="radio"/>    |
| <b>Units</b>                                   |  |                          |
| Use Regional Settings Units                    | (Enter only on Reference) Causes the New Study dialog box to set the patient Height and Weight measurement system according to regional settings (e.g., metric or U.S.). When not checked, the user can choose the desired system. | <input type="checkbox"/> |
| Link both Units                                | (Enter only on Reference) Only applicable when [Use Regional Settings Units] is not selected, this option forces the user to use the same measurement system for both patient height and weight (e.g., metric or U.S.).            | <input type="checkbox"/> |
| <b>Edit Patient Information</b>                |  |                          |
| Allow user to modify patient information       | (Enter only on Reference) If check this box, enable to modify study information without input administrator's password.  | <input type="radio"/>    |
| Exchange positions of Weight and Height fields | If check this box, show study information in order of "Height" and "Weight".   | <input type="radio"/>    |


## ■ Fonts and Colors Tab

Fonts and Colors options are configured as follows.



**NOTE**

Font & Colors is performed only on the Reference computer.



| Item                 | Description  | Mod                   |
|----------------------|--|-----------------------|
| <b>Color Options</b> |  |                       |
| Arrow                | Set Arrow color of Annotation.   | <input type="radio"/> |
| Annotation Text      | Set Text color of Annotation.  | <input type="radio"/> |
| Line                 | Set Line color of Annotation.  | <input type="radio"/> |
| Curve                | Set Curve color of Annotation.   | <input type="radio"/> |
| Caliper              | Set Caliper color of Annotation.   | <input type="radio"/> |
| Analysis Text        | Set Text color of QCA/LV.  | <input type="radio"/> |
| Study Text           | Set Text color of Study of Annotation.   | <input type="radio"/> |
| Change               | Click the [Change] button next to each item to open the Color window. Select a color for that item and click [OK]. | <input type="radio"/> |
| Default              | Click [Default] to reset color for all items to white.   | <input type="radio"/> |
| <b>Font Options</b>  |  |                       |
| Annotation Text      | Set font and font size of Annotation.  | <input type="radio"/> |
| Analysis Text        | Set font and font size of QCA/LV.  | <input type="radio"/> |


| Item                  | Description   | Mod                   |
|-----------------------|---|-----------------------|
| Study Text            | Set font and font size of Study of Annotation.  | <input type="radio"/> |
| Change                | Click the [Change] button next to each text type to open the Font window. Select a font, font style and size, and click [OK]. | <input type="radio"/> |
| Boxed                 | Select [Boxed] to have a box surrounding the specified type of text.  | <input type="radio"/> |
| Underline             | Select [Underline] to underline the specified type of text.   | <input type="radio"/> |
| Default               | Click [Default] to reset all text parameters to default.  | <input type="radio"/> |
| <b>Other Settings</b> |   |                       |
| Arrow Weight          | Defines the arrow weight.   | <input type="radio"/> |
| Caliper Step          | Defines the caliper step.   | <input type="radio"/> |
| Default               | Click [Default] to reset both arrow weight and caliper step parameters to default.  | <input type="radio"/> |
| <b>Image Selector</b> |   |                       |
| Font size             | Select [Font size] to set font size of angle information displayed on the image selector.                                     | <input type="radio"/> |
| Default               | Click [Default] to reset font size.   | <input type="radio"/> |


**NOTE**

When changing Font size of Image Selector on Acquisition monitor during study, the change will be available in the next study.




## Database Tab

The Database tab contains up to eight buttons, one for each table.




| Item   | Description   | Mod        |            |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |
|--|---|------------|------------|-----------|------|---|-------|-------|-------|--|----------|----------|------|--|-------|-------|-------|--|------------|------------|-------|--|----------|----------|-------|--|----------|----------|-----------|---|
| DICOM Tables   | The four DICOM tables are for DICOM experts only and are to be left as is.  | ×          |            |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |
| User Table   | <p>(Present only when Security feature is not installed.) This table enables the administrator to view, add, edit, and delete users. On the Database tab, click [User Table] to view all users. The table contains several items (columns) per user (row).</p> <div><p>User</p><table><thead><tr><th></th><th>User Login</th><th>User Name</th><th>Type</th></tr></thead><tbody><tr><td>▷</td><td>Admin</td><td>admin</td><td>SUPER</td></tr><tr><td></td><td>operator</td><td>operator</td><td>USER</td></tr><tr><td></td><td>super</td><td>super</td><td>SUPER</td></tr><tr><td></td><td>serviceapp</td><td>serviceapp</td><td>SUPER</td></tr><tr><td></td><td>shimadzu</td><td>shimadzu</td><td>SUPER</td></tr><tr><td></td><td>adminope</td><td>adminope</td><td>SUPERUSER</td></tr></tbody></table><div><div>Add/Modify</div><div>Delete</div><div>Update</div><div>Refresh</div><div>Close</div><div>Set Password</div></div></div> |            | User Login | User Name | Type | ▷ | Admin | admin | SUPER |  | operator | operator | USER |  | super | super | SUPER |  | serviceapp | serviceapp | SUPER |  | shimadzu | shimadzu | SUPER |  | adminope | adminope | SUPERUSER | × |
|  | User Login  | User Name  | Type       |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |
| ▷  | Admin   | admin      | SUPER      |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |
|  | operator  | operator   | USER       |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |
|  | super   | super      | SUPER      |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |
|  | serviceapp  | serviceapp | SUPER      |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |
|  | shimadzu  | shimadzu   | SUPER      |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |
|  | adminope  | adminope   | SUPERUSER  |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |
| <p>To add a user, click [Add]. A new blank row appears. Fill in each item in the row. To edit an item, click it and make the desired changes. To save changes, click a different row or [Add] or [Close]. The [Refresh] button also saves any changes and moves the pointer to the first row. To delete a user, select the user (click anywhere in the row) and click [Delete]. Do not delete standard users such as "super", "operator", and "serviceapp".</p> <p>See  <a href="#">"17.3 Viewing User Type" P.17-3</a></p> |   |            |            |           |      |   |       |       |       |  |          |          |      |  |       |       |       |  |            |            |       |  |          |          |       |  |          |          |           |   |




| Item   | Description   | Mod          |             |                  |             |           |        |      |                     |           |         |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
|--|---|--------------|-------------|------------------|-------------|-----------|--------|------|---------------------|-----------|---------|--------|--------|---|--|-----------|--|-----------|--|---|--------|------------|---------------------|-----------|--------|---|--|--------------|--|------------------|--|---|----|--------|----------|-----------|---------|---|--------|------------|-------------|-----------|--------|---|--|-------|--|-----------|--|---|
| Performing, Referring Physicians and Operator Tables | <p>(Reference monitor only) These tables enable the administrator to view, add, edit, and delete physicians. The physicians defined in these tables are made available for selection when creating a new study.</p> <p>See  <a href="#">"4.3 Entering Study Information" P.4-10</a>.</p> <p>On the Database tab, click [Performing Physician], [Referring Physician], or [Operator] to view physicians. The table contains separate fields (columns) for each part of the physician name with one physician per row.</p> <div><p><b>Performing</b></p><table><tr><th></th><th>Prefix</th><th>First Name</th><th>Middle Name</th><th>Last Name</th><th>Suffix</th></tr><tr><td>▷</td><td></td><td>Check</td><td></td><td>Daily</td><td></td></tr><tr><td></td><td></td><td>Second</td><td></td><td>Physician</td><td></td></tr></table><p><u>A</u>dd/<u>M</u>odify   <u>D</u>el<u>e</u>te   <u>U</u>p<u>d</u>ate   <u>R</u>ef<u>r</u>esh   <u>C</u>lose</p></div> <div><p><b>Referring</b></p><table><tr><th></th><th>Prefix</th><th>First Name</th><th>Middle Name</th><th>Last Name</th><th>Suffix</th></tr><tr><td>▷</td><td></td><td>Check</td><td></td><td>Daily</td><td></td></tr><tr><td></td><td></td><td>Second</td><td></td><td>Physician</td><td></td></tr></table><p><u>A</u>dd/<u>M</u>odify   <u>D</u>el<u>e</u>te   <u>U</u>p<u>d</u>ate   <u>R</u>ef<u>r</u>esh   <u>C</u>lose</p></div> <div><p><b>Operator</b></p><table><tr><th></th><th>Prefix</th><th>First Name</th><th>Middle Name</th><th>Last Name</th><th>Suffix</th></tr><tr><td>▷</td><td></td><td>First</td><td></td><td>Physician</td><td></td></tr></table><p><u>A</u>dd/<u>M</u>odify   <u>D</u>el<u>e</u>te   <u>U</u>p<u>d</u>ate   <u>R</u>ef<u>r</u>esh   <u>C</u>lose</p></div> <p>To add a physician, click [Add]. A new blank row appears. Fill in each part of the name. Prefix and Suffix can be left blank. At least enter the first and last names. To edit an item, click it and make the desired changes. To save changes, click a different row or [Add] or [Close]. The [Refresh] button also saves any changes and moves the pointer to the first row. To delete a physician, select the physician (click anywhere in the row) and click [Delete]. To modify a physician, delete the physician and add a new one.</p> |              | Prefix      | First Name       | Middle Name | Last Name | Suffix | ▷    |                     | Check     |         | Daily  |        |   |  | Second    |  | Physician |  |   | Prefix | First Name | Middle Name         | Last Name | Suffix | ▷ |  | Check        |  | Daily            |  |   |    | Second |          | Physician |         |   | Prefix | First Name | Middle Name | Last Name | Suffix | ▷ |  | First |  | Physician |  | ○ |
|  | Prefix  | First Name   | Middle Name | Last Name        | Suffix      |           |        |      |                     |           |         |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
| ▷  |   | Check        |             | Daily            |             |           |        |      |                     |           |         |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
|  |   | Second       |             | Physician        |             |           |        |      |                     |           |         |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
|  | Prefix  | First Name   | Middle Name | Last Name        | Suffix      |           |        |      |                     |           |         |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
| ▷  |   | Check        |             | Daily            |             |           |        |      |                     |           |         |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
|  |   | Second       |             | Physician        |             |           |        |      |                     |           |         |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
|  | Prefix  | First Name   | Middle Name | Last Name        | Suffix      |           |        |      |                     |           |         |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
| ▷  |   | First        |             | Physician        |             |           |        |      |                     |           |         |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
| Patient Table  | <p>(Reference monitor only) The Patient table enables the administrator to view and delete patients. The patients contained in this table are made available for selection when creating a new study without worklist management.</p> <p>See  <a href="#">"Define New Study with Local Worklist Management" P.4-15</a>.</p> <p>On the Database tab, click [Patient Table] to view patients. The table contains separate fields (columns) for each patient information item with one patient per row.</p> <div><p><b>Patient</b></p><table><tr><th></th><th>Prefix</th><th>First Name</th><th>Middle Name</th><th>Last Name</th><th>Suffix</th><th>Day</th><th>Month</th><th>Year</th><th>ID</th><th>Weight</th><th>Height</th></tr><tr><td>▷</td><td></td><td>Patient==</td><td></td><td>Name==</td><td></td><td>0</td><td>0</td><td>0</td><td>PICD:AP9500REV:0136</td><td>0</td><td>0</td></tr><tr><td></td><td></td><td>HANAKO=花子=11</td><td></td><td>SHIMADZU=島津=2721</td><td></td><td>8</td><td>12</td><td>1992</td><td>28504711</td><td>45.000000</td><td>165.000</td></tr></table><p>◀   <u>A</u>dd   <u>D</u>el<u>e</u>te   <u>R</u>ef<u>r</u>esh   <u>C</u>lose</p></div> <p>To delete a patient, select the patient (click anywhere in the row) and click [Delete]. The [Refresh] button moves the pointer to the first row. The Add feature of Patient Table is not enabled. Instead, add new patients when creating new studies in direct-entry mode.</p> <p>See  <a href="#">"Define New Study with Local Worklist Management" P.4-15</a>.</p>  |              | Prefix      | First Name       | Middle Name | Last Name | Suffix | Day  | Month               | Year      | ID      | Weight | Height | ▷ |  | Patient== |  | Name==    |  | 0 | 0      | 0          | PICD:AP9500REV:0136 | 0         | 0      |   |  | HANAKO=花子=11 |  | SHIMADZU=島津=2721 |  | 8 | 12 | 1992   | 28504711 | 45.000000 | 165.000 | ○ |        |            |             |           |        |   |  |       |  |           |  |   |
|  | Prefix  | First Name   | Middle Name | Last Name        | Suffix      | Day       | Month  | Year | ID                  | Weight    | Height  |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
| ▷  |   | Patient==    |             | Name==           |             | 0         | 0      | 0    | PICD:AP9500REV:0136 | 0         | 0       |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |
|  |   | HANAKO=花子=11 |             | SHIMADZU=島津=2721 |             | 8         | 12     | 1992 | 28504711            | 45.000000 | 165.000 |        |        |   |  |           |  |           |  |   |        |            |                     |           |        |   |  |              |  |                  |  |   |    |        |          |           |         |   |        |            |             |           |        |   |  |       |  |           |  |   |

## Devices Tab

Device options are configured as follows.




| Item                   | Description   | Mod |
|------------------------|---|-----|
| <b>Media Write</b>     |   |     |
| Write Speed            | Choose [Maximum] to write at the maximum speed supported by the writer and detected media. Choose any other value to reduce the maximum write speed. If you choose an unsupported speed, the speed is automatically adjusted. <b>You must restart the system after changing this value.</b>   | ×   |
| Convert to 512x512     | Defines whether studies are converted to a lossless compressed DICOM format (512x512 8-bit) before writing to removable media. If unchecked, the DICOM format (1k x 1k 12-bit) is written to the media. Large uncompressed Cine (raw) studies may not fit on a single CD so it is recommended to check this box if using CD.  | ○   |
| Automatically Write    | (Reference monitor only) When a study that is open for acquisition is closed, the user is automatically prompted to insert blank media onto which the study will be written. Check the media type for which the user should be prompted. If you check neither CD nor DVD, no automatic prompting occurs. If you check both, the user will be prompted to choose the media type. | ○   |
| Auto-accept Defaults   | Causes only the prompts related to inserting blank media to be displayed in a pop-up window. All other media writing-related prompts appear in the status bar of the Studies Management and Image Viewer windows. Furthermore, the prompt to select specific series is skipped.   | ○   |
| Image Processing       | Causes image processing including DSA subtraction to be applied to the images before writing them to CD/DVD. See  | ○   |
| Additional Files       | Enables you to include additional files (such as a player) on the media. Click[ Select] and choose a folder that contains the additional files. All files and folders in the selected folder will be written to the root of the media. File names must have the 8.3 format. Folder names must be no more than 8 characters in length.   | ×   |
| <b>Monitor Power</b>   |   |     |
| Disable Second Monitor | Not used.   | ×   |

| Item                         | Description  | Mod |
|------------------------------|--|-----|
| Turn Off                     | <p>When selected, the [After (Minutes)] value defines how many minutes (1 to 60) must elapse without keyboard, mouse, or IVR NEO/IVR Shuttle activity, before the monitor screen blanks and switches into its power-save mode. Windows screen savers must be disabled.</p> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;"> <div style="text-align: center;">  <b>NOTE</b> </div> <p>Turn Off feature is disabled whenever a study is open for acquisition or review.</p> </div> | ○   |
| Screen Saver                 | When selected, the [After (Minutes)] value defines how many minutes (1 to 60) must elapse without any activity, screen saver activates. When selecting [Password Prompt], you will be asked for a password when canceling screen saver.  | ○   |
| Mouse Wheel Scroll Direction | Set scroll direction of mouse wheel. When check [Invert Direction] and scroll the mouse wheel to downward to show the next frame.  | ×   |

External Software Tab


External Software options are configured as follows.



| Item           | Description   | Mod |
|----------------|---|-----|
| Shell Commands | Enter any shell command (with drive letter and path) that is to be executed before [Application Startup] or after [Application Exit]. For example, "c:\Shutdown.exe" could be entered in [Application Exit] to cause the system to shut down upon exit. | ×   |

## **Fusion Tab**

Use [Fusion] tab for setting communication with 3D workstation. Do not change the settings except Shimadzu or specified service personnel.



## **Logs Tab**


The [Log] tab enables you to log various events. Do not use this tab to activate logging unless you have been directed to do so by Shimadzu technical support.

## Cyber Security Tab



Cyber Security is initially disabled. Please contact our service representative to enable it.

Cyber Security options are configured as follows.




| Item                    | Description  | Mod                   |
|-------------------------|--|-----------------------|
| <b>User Account</b>     |  |                       |
| User Account Management | Create, edit and delete user account.  | <input type="radio"/> |
| Policies                | Set user ID and password policies.   | <input type="radio"/> |
| <b>Audit Log</b>        |  |                       |
| Audit Log Management    | Browse, export and delete audit log.   | <input type="radio"/> |
| Stored days             | Set the number of days to store audit log. (1 to 180)<br>Default is 90 (days). | <input type="radio"/> |

## User Account Management


Create, edit and delete user account.

Click [User Account Management] on the Cyber Security Configuration window. The User Account Management window appears.



- **Create User Account**

- 1 Click [Create].
- 2 The Create window appears.



### 3 Enter User ID and Password, and set User Level.

The screenshot shows a 'Create' dialog box with a dark background. It contains three input fields: 'User ID' with the value 'GuestUser', 'Password' with eight asterisks, and 'User Level' with a dropdown menu showing 'USER' selected. Below the dropdown menu, the options 'USER' and 'SUPER' are visible. At the bottom right, there are 'OK' and 'Cancel' buttons.

| Item                        | Description  |
|-----------------------------|--|
| USER                        | <p>Starting and ending a study</p> <p>Operation on the Studies Management window</p> <ul style="list-style-type: none"> <li>• Display of past study</li> <li>• Searching study</li> <li>• Browsing study information</li> <li>• Viewing RDSR</li> <li>• Sending study</li> <li>• Receiving study</li> <li>• Media writing of study information</li> <li>• Protection and unprotection of study</li> <li>• Deleting study</li> </ul>    |
| SUPER                       | <p>Functions that can be performed at the USER level</p> <p>Audit log related operations</p> <ul style="list-style-type: none"> <li>• Viewing audit logs</li> <li>• Outputting audit logs</li> </ul> <p>User Account Management Functions</p> <ul style="list-style-type: none"> <li>• Creation of user accounts</li> <li>• Deleting user accounts</li> <li>• Editing user accounts</li> <li>• User account policy settings</li> </ul> |
| (Available without Sign in) | <p>Operating side menu items</p> <ul style="list-style-type: none"> <li>• Fluoroscopy</li> <li>• [Funct.] menu</li> <li>• FPD Calibration</li> <li>• Error Report</li> <li>• Monitor</li> <li>• Configuration</li> <li>• Error List</li> <li>• Network</li> </ul>  |



- **Edit User Account**

- 1 Click [Edit].
- 2 The Edit window appears.

The 'Edit' window contains the following fields and controls:

- User ID:** A text box containing 'guestuser'.
- Password:** A section with two radio buttons: 'Keep existing password' (unselected) and 'Create new password' (selected). Below the radio buttons is an empty password input field.
- Checkboxes:** A checkbox labeled 'User must change password at next sign in' is checked.
- User Level:** A dropdown menu currently showing 'USER'.
- Buttons:** 'OK' and 'Cancel' buttons are located at the bottom right.

- 3 To change the password, select [Create new password] and register the new password.  
Also, if you check the "User must change password at next sign in" check box, the user can change the password at the first sign-in.

**NOTE**

Cannot change User ID.

- 4 Change User Level if necessary.

- **Delete User Account**

- 1 Select an User ID to delete, then click [Delete].

The 'Delete' dialog box contains the following text and controls:

- Title:** Delete
- Message:** Are you sure you want to delete the following user ID: <guestuser>?
- Buttons:** 'OK' and 'Cancel' buttons are at the bottom right.

- 2 A confirmation dialog appears and click [OK].  
User ID is deleted and the list will be updated.

• Policies

Set policies for User ID and Password.

Click [Policies] on the Cyber Security Configuration window. The Policies window appears.

**Policies**

**User ID**

Minimum Length (1-63)

**Password**

Minimum Length (0-63)

History Numbers (0-15)

Difference of password length from last password (0-15)

Type of password characters

Lower case letters ☐

Capital letters ☐

Numbers ☐

Special Characters ☐

Prohibited changing period (days) (0-15)

Expiration date (days) (0-365)

Password failure times for user account locking (0-20)

| Item   | Description  | Default Value<br>(Input Range) |
|--|--|--------------------------------|
| <b>User ID</b>                                   |  |                                |
| Minimum Length                                   | Set the minimum number of characters for user ID.  | 1<br>(1 to 63)                 |
| <b>Password</b>                                  |  |                                |
| Minimum Length                                   | Set the minimum number of characters for password.<br>If "0" is set, no password is required.  | 0<br>(0 to 63)                 |
| History Numbers                                  | Set the number of times password are recorded.<br>If "0" is set, setting the same password as the current password is permitted when changing the password.<br>If "1" is set, the same password as the current password is not permitted when changing the password. | 0<br>(0 to 15)                 |
| Difference of password length from last password | When changing password, set the difference of password length from last password.<br>If "0" is set, setting the same password as the current password is permitted.  | 0<br>(0 to 15)                 |

| Item  | Description   | Default Value<br>(Input Range)  |
|---|---|---|
| Type of password characters                     | Set the type of characters used for password.   | Uncheck<br>Capital letters (A to Z)<br>Lower case letters (a to z)<br>Numbers (0 to 9)<br>Special Characters (~ ! @ # \$ % ^ & * ( ) _ - + = { } [ ] \   : ; " < > , . ? /) |
| Prohibited changing period (days)               | Set the number of days that the password can be changed after changing the password.<br>If "0" is set, the password is allowed to be changed again even on the day it is changed. | 0<br>(0 to 15)  |
| Expiration date (days)                          | Set password expiration date.<br>If "0" is set, password will not expire.<br>If "90" is set, you will be requested for a new password on day 91.                                  | 0<br>(0 to 365)   |
| Password failure times for user account locking | Set the number of failed password entries that result in the user account is locked.<br>If "0" is set, the password is not locked.  | 0<br>(0 to 20)  |

### Audit Log Management

Enable to browse, export, and delete audit log on audit log management.




Deleting audit log is only available in serviceapp.

Leave the following items in the audit log:


- Sign in / Sign out
- Operation of Studies Management window
  - Show past study
  - Search for study
  - Browse study
  - Viewing RDSRs
  - Sending study
  - Receiving study
  - Outputting study to Media
  - Protecting study / Unprotecting study
  - Deleting study
  - Anonymizing study
- Operations related to audit logs
  - Viewing audit log
  - Output of audit log
  - Delete of audit log
- User account management functions
  - Creating user accounts
  - Editing user accounts
  - Deleting user accounts
  - Editing user account policies

- **Browse Audit Log**

- 1 Click [Audit Log Management] on the Cyber Security Configuration window.
- 2 The Audit Logs window appears.




- 3 Select a date and audit logs are displayed.



- **Export Audit Log**

- 1 Select a date to export. Multiple selections are available.
- 2 Click [Export].



- 3 The Media Writing window appears. Click [Start].



- 4 Start writing media.

### 17.7.4 Menus and DUP Configuration


**NOTE**

Any changes you make to menus and DUPs do not take effect until you create a new study.

**NOTE**

Menus and DUP configuration is performed only on the Reference computer.

The [Menu Creation] group of the [Menus Configuration] tab enables the installation personnel and administrator to configure Fluoro, Rad, and Comment menus, and Fluoro and Rad Digital User programs (DUPs).



In the [Menu Category] list you can choose [Comments], [Fluoro], or [Rad] to configure the respective menus.

## ■ Configure the Comment Menu


In DAR-9500f, the Comment menu can be displayed whenever a study is open for acquisition by clicking [Comm] on the side menu. The menu displays a list of standard comments about various aspects of the procedures performed.

The content of this menu is site-specific and is typically initially configured by the installation personnel according to customer wishes.

To configure the Comment menu, follow this procedure.

- 1 In the [Menu Category] list, choose [Comments].

The current Comments menu configuration is show below.




- 2 To add a new comment (up to 50 total), click [New].  
The comment "New Comment" appears at the bottom of the list.
- 3 Select the [New Comment] item in the list, and then in the [Current Selection] box change the text "New Comment" to whatever you like.  
(1 to 48 characters, excluding underscore "\_", back slash "\", number "#", equals "=", angle brackets "<>", and ":" colon)
- 4 Click [Update].  
The comment list is updated with the new text of your comment.
- 5 Position the comment with [Up/Down] as desired.
- 6 To sort comments in alphabetical order, click the [Sort] button.
- 7 To delete a comment, select it and then click [Delete].
- 8 When finished, click [Apply] to save your changes.



■ **Configure Fluoro/Rad Menus and DUPs**

In DAR-9500f, Fluoro and Rad menus enable the operator to quickly choose digital user programs (DUPs) that define exposure and acquisition characteristics. For example, a Rad menu, with its two levels shown.



The Menu creator enables you to define menu item text and sequence for both menu levels, plus DUP parameters for second-level menu items. The menus created here are displayed by the operator when preparing to perform an acquisition. The content of these menus and associated DUP parameters are site-specific and are typically configured by the installation personnel according to customer wishes.


■ **Configure Fluoro/Rad Menus**

To configure Fluoro and Rad menus follow this procedure (except where otherwise noted, this procedure applies equally to Fluoro and Rad).


- 1 In the [Menu Category] list, choose either [Fluoro] or [Rad].  
The current menu configuration (first level) is shown.



- 2 To see the second menu level, select a menu item and click [Edit].  
The second level menu appears.



- 3 Each item on second-level menus has a set of Digital User Program (DUP) parameters associated with it.  
For details on how to configure DUP parameters, see the next section "[Configure Fluoro/Rad DUPs](#)".
- 4 To position a menu item vertically, select it and then move it up or down with the [Up/Down] buttons.
- 5 To add a new menu item, first select an existing menu item whose DUP parameters most closely match the purpose of the new DUP and click [New].  
A copy of the item is added to the menu item list. Select the new (copied) item and edit its text in the Current Selection box (1 to 48 characters, excluding underscore "\_", back slash "\", number "#", equals "=", angle brackets "<>", and colon ":").
- 6 Click [Update]. Position the new item with [Up/Down] as desired.  
A limited number of menu items are permitted at each level. If you attempt to add too many, a warning message similar to this is displayed. Either delete an unnecessary menu item and then click [New] again, or rename / modify an existing item.



- 7 To rename a menu item, select it, enter the new name in the Current box and then click [Update].  
Although each menu can be up to 48 characters long, some characters will be truncated on the right side if the full menu name does not fit in one line. Therefore, after configuring menu names, verify that each name displays without truncation.
- 8 To delete a menu item, select it and then click [Delete].  
The first menu item at each level is always present and cannot be deleted, however, it can be renamed.
- 9 To go back to a previous menu level, click [Back].
- 10 When finished, click [Apply] to save your changes.


## Configure Fluoro/Rad DUPs



It is a good practice to work on a copy of an existing good DUP rather than directly modifying the DUP. To do this, select the DUP to be copied, click [New], adjust or enter a new name in Current Selection, and then click [Update]. Then edit the copy with your new values.

To configure Fluoro and Rad digital user programs (DUPs), follow this procedure (except where otherwise noted, this procedure applies equally to Fluoro and Rad).

- 1 Navigate (arrow keys or mouse) to the desired second-level menu item and click the [Edit] button.



The DUP parameters are displayed like this (Rad and then Fluoro shown):

<RAD Type: Other then HQ-DSA>

**RAD: Cardio/CAG[15f-10s]**

Program Description:  ☒ Frontal ☒ Lateral ☒ Biplane ☐ Timer

Camera

Frame Rate (fps):  Rad Time:

☒ Double Speed Acquisition

Grab Resolution:  Rad Injector Synch:

Exposure Time (ms):  Rad I-A (s):

Rad Type:  Rad M-A (s):

Rad Mask:  Body Part Examined:

IOF Level:  STENTVIEW

FPD Mode No:

Application:

App. Param:

Generator

Condition Number:  BH Filter:  Density:

Max. Number of Images:  Rad Mode:

Image Pre-Processing

☒ EnhancedView

Pre-Processing:

LUT:

Image Post-Processing

☒ DSA ☐ 3DWS ☐ 3D Recon Mode:

Normal ☒ DSA

Edge Enhancement:  Power:  Reverse Video:

AWL:  Brightness:  Contrast:

OK Cancel Apply

<RAD Type: HQ-DSA>

**RAD: Head/DSA[30s]**

Program Description:  ☒ Frontal ☒ Lateral ☒ Biplane ☐ Timer

Camera

State Acquisition:      Total:

Rad Time:

Grab Resolution:  Rad Injector Synch:

Exposure Time (ms):  Rad I-A (s):

Rad Type:  Rad M-A (s):

Rad Mask:  Body Part Examined:

IOF Level:  STENTVIEW

FPD Mode No:

Application:

Generator

Condition Number:  BH Filter:  Density:

Max. Number of Images:  Rad Mode:

Image Pre-Processing

☒ EnhancedView

Pre-Processing:  Temporal DSA:  Power:

LUT:

Image Post-Processing

☒ DSA ☐ 3DWS ☐ 3D Recon Mode:

Normal ☒ DSA

Edge Enhancement:  Power:  Reverse Video:

AWL:  Brightness:  Contrast:

OK Cancel Apply

The screenshot shows the 'FLUORO: Head/10pps' window with the following settings:

- Program Description:** Fluo
- Camera:**
  - Frame Rate (fps): 10
  - Grab Resolution: Grab FPD Fluoro 1024 12bit
  - Exposure Time (ms): 16
  - Fluoro RSM: RSM-0
  - Record Rate (fps): 0
  - Max Record Time (s): 20
  - Fluoro Dose: Normal
  - IOF Level: Normal
  - FPD Mode No: 8
  - Body Part Examined: HEAD
- Generator:**
  - Condition Number: 13
  - BH Filter: 0
  - Density: 0
- Image Pre-Processing:**
  - Temporal: K=1
  - Temporal DSA: K=1\_125
  - LUT: None
  - Parameter:
    - BG Control: BG 1
    - LUT: Vessel 2, Power: 120
    - DSA-MAP: Vessel 2, Power: 120
    - FluoroMAP: Vessel 2, Power: 120
- Image Post-Processing:**
  - DSA: Normal
  - Edge Enhancement: Convo-06
  - AWL: AWL 07
  - Power: 0
  - Brightness: 24
  - Contrast: -4
  - Reverse Video: On

Buttons at the bottom: OK, Cancel, Apply.

## 2 Edit the DUP parameters as needed.

The purpose of each parameter is briefly described in the following table (differences between Fluoro and Rad parameters are indicated).



DUP parameters only set the initial values to be used during subsequent acquisitions.

| No            | Group/Item               | Description  | R                     | F                     | Mod                   |
|---------------|--------------------------|--|-----------------------|-----------------------|-----------------------|
| 1             | Program Description      | Defines the Acquisition Active message that is displayed in the upper-right area of the Image Viewer window during actual image acquisition (e.g., "DA", "RSM", or "Fluo"). During Direct Fluoro Record, the additional word "rec" is added to the right as in "Fluo rec". (Any characters, except underscore "_", back slash "\", number "#", equals "=", and angle brackets "<>"). | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 2             | Timer                    | The timer will start when acquisition starts. Use Timer button on the side menu or IVR NEO/IVR Shuttle/SMART Touch button to stop the timer.   | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| <b>Camera</b> |                          |  |                       |                       |                       |
| 3             | Frame Rate               | The X-ray pulse rate in frames per second.   | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 4             | Double Speed Acquisition | Not used.  | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| 5             | Grab Resolution          | The acquisition image resolution (1024x1024 or 512x512) and bit depth.   | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |

| No        | Group/Item                 |                          | Description   | R                     | F                     | Mod                   |
|-----------|----------------------------|--------------------------|---|-----------------------|-----------------------|-----------------------|
| 6         | Exposure Time              |                          | The exposure time in milliseconds for the acquisition of one image.   | <input type="radio"/> | <input type="radio"/> | ×                     |
| 7         | Acquisition Plane (F/L/Bi) |                          | The plane for which acquisition will occur: Frontal, Lateral, or Both   | <input type="radio"/> |                       | ×                     |
| 8         | Rad Time                   |                          | The total radiographic time in seconds, not including Spot.   | <input type="radio"/> |                       | ×                     |
| 9         | Rad Type                   |                          | The radiography type such as DA or DSA.   | <input type="radio"/> |                       | ×                     |
| 10        | Rad Mask                   |                          | For FPD Rad Mask. Sets the number of frames in the mask. (1 to 16)  | <input type="radio"/> |                       | ×                     |
| 11        | Rad Injector Synch         |                          | Turns injection synchronization ON or OFF.  | <input type="radio"/> |                       | ×                     |
| 12        | Rad I-A                    |                          | The delay in seconds between injection and beginning of live acquisition.   | <input type="radio"/> |                       | <input type="radio"/> |
| 13        | Rad M-A                    |                          | The delay in seconds between mask acquisition and beginning of live acquisition.  | <input type="radio"/> |                       | <input type="radio"/> |
| 14        | Fluoro RSM                 |                          | RSM filter type: [RSM-0], [RSM-1], [RSM-2]  |                       | <input type="radio"/> | ×                     |
| 15        | Record Rate                |                          | The maximum rate (frames per second) at which Fluoro can be recorded. (0 to Frame Rate (30 for Continuous).)  |                       | <input type="radio"/> | ×                     |
| 16        | Max Record Time            |                          | The maximum number of seconds that Fluoro can be recorded for one acquisition.  |                       | <input type="radio"/> | <input type="radio"/> |
| 17        | Fluoro Dose                |                          | Dose mode: Normal, Low, or High.  |                       | <input type="radio"/> | ×                     |
| 18-1      | Stage Acquisition          | Exposure Time            | Rad Time (sec) for each frame rate. If "0" is set, the frame rate will not be used.   | <input type="radio"/> |                       | <input type="radio"/> |
| 18-2      |                            | Total Number of Exposure | Total number of exposure of all frame rate. For 3fps and over acquisition rate, total number of exposure is limited to 90. For all frame rate, total number of exposure is limited to 150. No exposure will be performed for more than those limited number of times. | <input type="radio"/> |                       | <input type="radio"/> |
| 19        | IOF Level                  |                          | Not used.   | <input type="radio"/> | <input type="radio"/> | ×                     |
| 20        | FPD Mode No.               |                          | Specifies the acquisition table number. Usually do not modify.  | <input type="radio"/> |                       | ×                     |
| 21        | Sub                        |                          | Not used.   | <input type="radio"/> |                       | ×                     |
| 22        | Body Part                  |                          | Selects the body part where to be studied.  | <input type="radio"/> |                       | ×                     |
| 23        | Application                |                          | Set an application (option) to be used.   | <input type="radio"/> |                       | ×                     |
| Generator |                            |                          |   |                       |                       |                       |
| 24        | Condition number           |                          | Generator Program Number. (0 to 255)  | <input type="radio"/> | <input type="radio"/> | ×                     |
| 25        | Max. number of Images      |                          | Specifies the maximum number of images (frames) that can be acquired during one acquisition for this program, regardless of how long the exposure switch is pressed. (0 to 1023)  | <input type="radio"/> |                       | <input type="radio"/> |
| 26        | BH Filter                  |                          | Specify the low energy X-ray reduction filter. "0" is the system standard.  | <input type="radio"/> |                       | ×                     |
| 27        | Rad Mode                   |                          | Not used.   | <input type="radio"/> |                       | ×                     |


| No                    | Group/Item               |              | Description   | R                     | F                     | Mod |
|-----------------------|--------------------------|--------------|---|-----------------------|-----------------------|-----|
| 28                    | Density                  |              | Specify Input Dose for FPD.<br>Change approximately 15% in one step.<br>“0” is the system standard.   | <input type="radio"/> |                       | ×   |
| Image Pre-Processing  |                          |              |   |                       |                       |     |
| 29                    | EnhancedView             |              | This function is not available.   | ×                     | ×                     | ×   |
| 30                    | Pre-Processing           |              | Enables image pre-processing filters such as Temporal, Temporal DSA and Fluoro Map. Choose [None] to disable pre-processing. Image pre-processing is applied to the images as they are being acquired, before saving.   | <input type="radio"/> | <input type="radio"/> | ×   |
| 31                    | Filter                   | Temporal     | Temporal filtering can be selected to reduce image noise with a sliding-window-averaging technique. The [Power] value (-100 to 100) controls filter behavior.<br>A value of -100 means only the previous temporal result is used, 0 means 50 % of the previous temporal result plus 50 % of the current image are used, 100 means only the current image is used (same as no temporal filtering). Useful range is from -20 to 90.<br>Each time you choose an item from the list, the corresponding default [Power] value is selected. | <input type="radio"/> | <input type="radio"/> | ×   |
| 32                    |                          | Temporal DSA | As same as [Temporal], but optimize for FluoroMAP.  |                       | <input type="radio"/> | ×   |
| 33                    | LUT                      |              | The [LUT] list provides several predefined LUT values used to control how images are compensated to achieve better quality. [Reverse Video] can be enabled if desired. Each time you choose an item from the [LUT] list, the corresponding default [Reverse Video] setting is selected.   | <input type="radio"/> | <input type="radio"/> | ×   |
| 34                    | BG Control               |              | Select parameters for the background image used for MAP function.   |                       | <input type="radio"/> | ×   |
| 35                    | Vessel LUT               |              | Select from the list setting LUT and a power apply to the mask image used for Fluoro MAP, SIMAP (LIVE) and SIMAP (Sub).   |                       | <input type="radio"/> | ×   |
| Image Post-Processing |                          |              |   |                       |                       |     |
| 36                    | Image Post-Processing    |              | Enables image post-processing that is applied to the image as they are being displayed. This processing generally does not change the image files.  | <input type="radio"/> | <input type="radio"/> | ×   |
| 37                    | DSA                      |              | Enables DSA image processing.   | <input type="radio"/> | <input type="radio"/> | ×   |
| 38                    | 3DWS                     |              | Enable the necessary image post-processing to prepare images for transmission to a 3D workstation. This post-processing makes image files that are compatible with the 3D workstation.  | <input type="radio"/> |                       | ×   |
| 39                    | 3D Recon Mode            |              |   | <input type="radio"/> |                       | ×   |
| 40                    | Filter, Edge Enhancement |              | Images can be made to appear sharper by applying Unsharp mask filtering or Convolution filtering. The Edge Enhancement list contains several predefined filters such as [Unsharp 1]. The [Power] value (-100 to 100) defines the strength of the selected filter. Negative values smooth the image, and positive values sharpens the image. Each time you choose an item from the [Edge Enhancement] list, the corresponding default [Power] value is selected.   | <input type="radio"/> | <input type="radio"/> | ×   |

| No | Group/Item            |                  | Description  | R                     | F                     | Mod |
|----|-----------------------|------------------|--|-----------------------|-----------------------|-----|
| 41 | Auto Window Level     | AWL              | Auto Window Level provides advanced automatic Brightness/Contrast optimization. The [AWL] list provides several predefined Auto Window Level configurations.<br><br>The AWL None choice corresponds to Auto Window Level not active (as when the Auto W. Level button on the Reference monitor side menu is not pushed.)<br><br>All other choices correspond to Auto Window Level Active (as when the [Auto W. Level] button on the Reference side menu is pushed.)<br><br>Each time you choose an item from the AWL list, the corresponding initial Brightness and Contrast values are selected.                            | <input type="radio"/> | <input type="radio"/> | ×   |
| 42 |                       | Brightness       | Set a number in the -100 to 100 range.   | <input type="radio"/> | <input type="radio"/> | ×   |
| 43 |                       | Contrast         | Set a number in the -100 to 2000 range (notice the different maximum than that of Brightness).   | <input type="radio"/> | <input type="radio"/> | ×   |
| 44 | Auto Window Level DSA | DSA AWL          | Auto Window Level DSA provides advanced automatic Brightness/Contrast optimization. The [AWL] list provides several predefined Auto Window Level DSA configurations. The [AWL None] choice corresponds to Auto Window Level DSA not active (as when the [Auto W. Level] button on the Reference monitor side menu is not pushed.)<br><br>All other choices correspond to Auto Window Level DSA Active (as when the [Auto W. Level] button on the Reference side menu is pushed.)<br><br>Each time you choose an item from the [AWL] list, the corresponding initial [DSA Brightness] and [DSA Contrast] values are selected. | <input type="radio"/> |                       | ×   |
| 45 |                       | DSA Brightness   | Set a number in the -100 to 100 range.   | <input type="radio"/> |                       | ×   |
| 46 |                       | DSA Contrast     | Set a number in the -100 to 2000 range (notice the different maximum than that of Brightness).   | <input type="radio"/> |                       | ×   |
| 47 |                       | Mask Weight      | Set the default of mask factor in subtraction.   | <input type="radio"/> |                       | ×   |
| 48 |                       | Mask Integration | Set the number of frames for integrating mask.   | <input type="radio"/> |                       | ×   |

R= Parameters apply to Rad DUP, F= Parameters apply to Fluoro DUP

### 3 Most of DUP parameters have default values.

Right-click DUP parameter to display [Reset to default] button, and click to reset to default value.



This page is intentionally left blank.



# *Chapter 18*

---

## Troubleshooting





This chapter provides information on what to do if trouble is encountered while using the system.

### Description

|                               |       |
|-------------------------------|-------|
| 18.1 Error Messages .....     | 18-2  |
| 18.2 Startup the System ..... | 18-5  |
| 18.3 Error Report .....       | 18-7  |
| 18.4 Unable to Sign In .....  | 18-10 |
| 18.5 Other .....              | 18-11 |

# 18.1 Error Messages

If you see one of the following error messages on the system monitor, take the suggested action to resolve the problem.

| No. | Message   | Action  |
|-----|---|---|
| 1   | An error occurred during cache space reservation.                                     | If you operate with CD storage, write any unsaved data to the CD. If not, please contact our service representative.  |
| 2   | No available cache.   |   |
| 3   | Burn Operation Failed.  | Insert a blank CD-R and re-execute the operation.   |
| 4   | Can not retrieve media information.   | Confirm whether the media is inserted properly.   |
| 5   | Cannot write DICOMDIR to CD-R (DVD-R). No patient has been added.                     | Confirm that writable media is inserted.  |
| 6   | Character <x> is not valid.<br>"<x>" represents variable information.                 | Try other characters.   |
| 7   | Could not find the PPS Manager. Please configure your PPS Manager in the hosts table. | Please contact our service representative.  |
| 8   | Error: Invalid query criteria entered.  | Input an appropriate value.   |
| 9   | No DICOM media found in drive <x>.  | Insert DICOM media.   |
| 10  | Cannot send MPPS of patient xxxx to MPPS manager xxxx.                                | <p>Check communication with MPPS manager. To cancel retry of MPPS support, click  at the bottom menu.</p>  <p>Select an error of MPPS support in the Error List and click [Manage Error]. And click [Yes] on a dialog.</p>   |



## Error Messages Related to FPD




When “FPD ERROR” is displayed at upper left of the system monitor, errors in the following list are also displayed. If an error message related to FPD is displayed, stop study immediately and contact Shimadzu service representative.

| No. | Message                           | Description  | Action  |
|-----|-----------------------------------|--|---|
| 1   | [Cooling Unit]: Water flow error  | Error of FPD cooling unit water flow.  | Stop study immediately and contact Shimadzu service representative. |
| 2   | [Sensor]: Temperature error       | Error of sensor chassis temperature.   | Stop study immediately and contact Shimadzu service representative. |
| 3   | [HV/DC Power]: Emergency shutdown | FPD shut down immediately. Fluoroscopy and radiography become disabled.                | Stop study immediately and contact Shimadzu service representative. |
| 4   | [Cooling FAN-1]: Error detected   | Error of FAN cooling ISYS and IOF board.   | Please contact our service representative.                          |
| 5   | [Cooling FAN-2]: Error detected   | Error of FAN cooling ISYS and IOF board from the back of a front panel of DAR cabinet. | Please contact our service representative.                          |





## Error Messages Related to SMART Touch


If an error message related to SMART Touch is displayed on the SMART Touch, contact Shimadzu service representative.

| No. | Message  | Action  |
|-----|--|---|
| 1   |  <p>The dialog box shows a red warning triangle icon. The text reads: "Connection Error. Connection with DAR-9500f digital system could not be established. This dialog will be hidden if reconnection is successful. Please press Restart if this message remains visible." There is a green "Restart" button with a circular arrow icon.</p> | If the Connection Error message is still displayed on the SMART Touch after reboot the application, reboot the system. If it does not solve the problem, contact Shimadzu service representative. |
| 2   |  <p>The dialog box shows a red warning triangle icon. The text reads: "Joystick Error. Joystick is not connected. Joystick is currently unavailable. SMART Touch functions other than joystick can be used." There is a grey "Close" button with an 'X' icon.</p>  | Touch [Close] and the SMART Touch functions other than joystick can be used. To use joystick, contact Shimadzu service representative.  |

| No. | Message  | Action   |
|-----|--|--|
| 3   |  <p>The dialog box titled "Joystick Error" features a red warning triangle icon on the left. The text reads: "Joystick is disconnected. Joystick is currently unavailable. SMART Touch functions other than joystick can be used." A "Close" button with a red 'X' icon is located in the bottom right corner.</p>  | Touch [Close] and the SMART Touch functions other than joystick can be used. To use joystick, contact Shimadzu service representative.                         |
| 4   |  <p>The dialog box titled "Joystick Error" features a red warning triangle icon on the left. The text reads: "An abnormal value was detected from the joystick. Please press Retry without touching the joystick. Pressing Close will closes this dialog. SMART Touch functions other than joystick can still be used." There are two buttons in the bottom right: "Retry" with a circular arrow icon and "Close" with a red 'X' icon.</p><br> <p>The dialog box titled "Joystick Validation" features a red warning triangle icon on the left. The text reads: "Please refrain from using the joystick while this message is visible."</p> | Touch [Close] and the SMART Touch functions other than joystick can be used. If error occurs though touching [Retry], contact Shimadzu service representative. |

# 18.2 Startup the System

| Problem   | Action   |
|---|--|
| <p>The power of monitor in the control room is not turned on, even turn ON the power of the system.</p>   | <p>Press the power of digital reference cabinet.</p>   |
| <p>A progress bar in [Configure Display Unit] dialog displayed on the Acquisition monitor does not proceed.</p>  | <p>Turn on the power supply of digital acquisition cabinet according to the following procedure.</p> <ol style="list-style-type: none"> <li>1 Press [Cancel] on the [Configure Display Unit] dialog.</li> </ol>  <ol style="list-style-type: none"> <li>2 Turn OFF the breaker of digital acquisition cabinet in the equipment room. And then wait approximately 60 seconds.</li> </ol>  |

| Problem | Action  |
|---------|---|
|         | <p><b>3</b> Turn ON the breaker of digital acquisition cabinet in the equipment room. Check that Cine Area Zoom operates normally.</p>  <p>The diagram shows a tall, white digital acquisition cabinet. A callout bubble points to a circuit breaker panel on the front. Inside the callout, a red circle highlights a switch, and a red arrow points upwards, indicating it should be turned ON. The panel also includes a blue indicator light and a black knob.</p> <ul style="list-style-type: none"> <li>• Cine Area Zoom</li> <li>• LIDM</li> </ul> |


# 18.3 Error Report



Making a contract with remote maintenance is necessary to enable Error Report.

Errors and screen capture will be reported when error occurs while using the system.

- 1 Click [Funct.]-[Report] on the side menu of REF monitor.  
[Error Report] window is displayed.



**Error Report**

1 — Q1. Please tell us about the error.  
Others



2 — Q2. Please tell us the current status.  
☐ Restored Unknown  
☐ Continues  
☒ Unknown

3 — Comments:

4 — ☒ Attach the current screen capture: Nb. Char 0 (<1024)

5 —   6


| No. | Item                                | Description   |
|-----|-------------------------------------|---|
| 1   | Q1. Please tell us about the error. | <div>Select an error from the followings:</div> <ul style="list-style-type: none"><li>Acquired image is abnormal</li><li>Troubles occurred on the C-Arm or the catheterization table</li><li>No X-rays are exposed</li><li>Others</li></ul> |

| No. | Item                               | Description  |
|-----|------------------------------------|--|
| 2   | Please tell us the current status. | <p>Select the current status from the followings:</p> <ul style="list-style-type: none"> <li>• Restored</li> <li>• Continuous</li> <li>• Unknown</li> </ul> <p>And if [Restored] is selected, select recovery procedure from the followings:</p> <ul style="list-style-type: none"> <li>• Cabinet breaker ON/OFF</li> <li>• Reboot the system</li> <li>• Reboot the application</li> <li>• Restore spontaneously</li> <li>• Unknown</li> </ul> |
| 3   | Comments:                          | Input error description.   |
| 4   | Attach the current screen capture: | <p>Check this item to capture the image displayed on ACQ and REF monitor. Include an annotation for screen capture.</p> <div style="border: 1px solid black; padding: 10px; margin: 10px 0;"> <div style="text-align: center;">  <b>NOTE</b> </div> <p>Captured screen may include personal data.</p> </div>  |
| 5   | Submit                             | Send an error report.  |
| 6   | Cancel                             | <p>Cancel input description and close [Error Report].</p> <div style="text-align: center;">  </div>  |

2 Input necessary information and click [Submit].



- 3 Confirmation message is displayed. Click [Yes].



Click [No] to cancel the sending of an error report.



An error report could not send if a message of [The connection with the server is not established yet.] is displayed after sending the report. Please contact our service representative.





For an emergency contact for serious trouble, please contact our service representative by phone.





An error report may not send due to the communication line trouble. If no response from our service representative, please contact by phone.

## 18.4 Unable to Sign In

If you are unable to sign in while the cybersecurity is enabled, do the following:

- Able to contact with the system administrator  
Have the administrator change the password.  
 ["• Edit User Account" P.17-63](#)
- Unable to contact with the system administrator  
 ["3.2.6 Emergency Sign In" P.3-8](#)

# 18.5 Other

| Problem  | Action  |
|--|---|
| The mouse pointer is locked within one mirror.   | Press the keyboard [Scroll Lock] key to de-activate the function.   |
| Unknown errors are displayed.<br> | Application is not responding.<br>Click [Yes] and restart the application.  |
| The screen became black.<br>Something like a grayscale image is displayed on the monitor.                          | Press the second button from the left on the lower-right corner of the monitor.<br>                |
| Unable to use SMART Touch.<br>Unable to use a joy stick of SMART Touch.  | Press the power supply button on the back of SMART Touch and turn off the power. And when the screen display of SMART Touch disappeared, turn on the power of SMART Touch to startup. |

This page is intentionally left blank.

# *Chapter 19*

---

## Maintenance

This chapter describes how to maintain the system and includes the following.

### Description

|                                 |      |
|---------------------------------|------|
| 19.1 Introduction .....         | 19-2 |
| 19.2 Daily Maintenance .....    | 19-3 |
| 19.3 Periodic Maintenance ..... | 19-5 |

## 19.1 Introduction

Maintenance of the system allows for vivid diagnostic images while assuring safety of the operator and patient.

The system is shipped in optimum condition through quality management and inspections. Regular maintenance is necessary to maintain the condition of the system.

The preventive inspection is classified as follows.

| No. | Type                 |                            | Performance            | Maintenance Cycle        |
|-----|----------------------|----------------------------|------------------------|--------------------------|
| 1   | Daily Maintenance    | Startup Maintenance        | Operator               | Daily (before operation) |
|     |                      | Post-operation Maintenance | Operator               | Daily (after operation)  |
| 2   | Periodic Maintenance | Calibration                | Operator               | Every 3 months           |
|     |                      | Preventive Maintenance     | Service Representative | Every 6 months           |

## 19.2 Daily Maintenance

The daily maintenance involves the startup and post-operation maintenance. The daily maintenance allows for longer duration of the product.

Startup maintenance procedures and finishing operations described in this section should be followed.

### 19.2.1 Startup Maintenance

#### **Purpose**

The startup inspection, prior to operations, is necessary for perfect and trouble-free performance during examination of patients.

#### **Procedures**

Follow the instruction below.

- 1 Visually check for the following.

| No. | Check Point                         | Procedure  |
|-----|-------------------------------------|--|
| 1   | Pinched, twisted or stripped cables | Fix or repair as required.   |
| 2   | Damage on chassis or cable          | Check visually and if something is wrong, contact Shimadzu Service Representative. |

- 2 Check the FPD Cooling System as follows. The FPD Cooling System is located in equipment room.

| No. | Check Point   | Procedure   |
|-----|---|---|
| 1   | Temperature displayed in front of the Cooling System must be $30 \pm 1^{\circ}\text{C}$ .   | If the displayed temperature is not within specified limits, contact Shimadzu Service Representative.   |
| 2   | There is enough water in the tank.  | Replenish with specified distilled water about once a month. (P/N:502-40086, Container for Distilled Water, 1L).<br>The appropriate solution level must be between the two level indicators.<br>The system does not operate if the solution level is below the low-level indicator. |
| 3   | Temperature and humidity of the examination and equipment room must be within the following:<br>Temperature:<br>10 to $35^{\circ}\text{C}$ (examination room)<br>10 to $35^{\circ}\text{C}$ (equipment room)<br>Humidity: 15 to 75% | If temperature or humidity is not within specified, condensation may occur. Adjust temperature of air-conditioning.   |



HEC002-A5B

- 3 Turn the power ON.
- 4 Verify that the system is operating normally by checking for the following.

| No. | Check Point   | Procedure   |
|-----|---|---|
| 1   | No error messages on both the Acquisition and Reference monitors. | If there is an error message, refer to Chapter 13 Troubleshooting.<br>If you are unable to resolve the problem, contact technical support.  |
| 2   | Fluoroscopy images are displayed.                                 | Make sure that fluoroscopy images are properly displayed on monitors.<br>If there is an error message, refer to Chapter 13 Troubleshooting.<br>If you are unable to resolve the problem, contact technical support.   |
| 3   | Radiography images are displayed.                                 | If a dummy patient does not exist, create one. Use the dummy patient to make sure that radiography images are properly displayed on monitors.<br>If there is an error message, refer to Chapter 13 Troubleshooting.<br>If you are unable to resolve the problem, contact technical support. |

## 19.2.2 Finishing Operations

This section describes the post-operation maintenance.

### ■ Post-operation Maintenance

Check the following when finishing operations.

| No. | Check Point                  | Procedure |
|-----|------------------------------|-----------|
| 1   | Dust or dirt is not removed. | Remove    |



## 19.3 Periodic Maintenance


Periodic maintenance includes calibration which is performed by a technician and preventive maintenance performed by the Shimadzu Service Representative.

### 19.3.1 Calibration

Automatic FPD calibration is supported to the system. Manual calibration is not necessary, so follow the instruction below when manually calibrating the system.

#### ■ FPD Calibration

- 1 Click [Funct]-[FPD] on the side menu.  
The Flat Panel Detector Calibration dialog box appears.
- 2 Make sure that the [Timeout] value is set to 10 seconds.  
This defines the maximum number of seconds to wait before canceling if there is a problem communicating with the FPD.
- 3 Select [Offset-All Modes] in the [Calibration Type] list.



- 4 When ready, click [Calibrate] button.  
Calibration progress messages are displayed. The Acquisition Mode item indicates the mode being calibrated. Completion messages are displayed when calibration is completed. Completion messages are displayed when calibration is completed. After the calibration, select the fluoroscopy mode which will be used for the first study.
- 5 If you must interrupt the calibration before it finishes, click [Abort] button, then click the [Close] button to close the window.  
You must later re-do the calibration before performing acquisitions.



Calibration data is not updated if you click the [Cancel] button.

### 19.3.2 Preventive Maintenance

Preventive maintenance is performed by the Shimadzu Service Representative. This maintenance requires sufficient knowledge of the internal mechanism and is hazardous.

The maintenance is performed every 6 months and will be charged after the guarantee expires.

#### ■ Periodic Part Replacement

To maintain system performance, some parts must be periodically replaced. Part replacement must be performed by a Shimadzu Service Representative.

Parts that require periodic replacement are shown as follows:

| Part                                       | Part No.     | Rated Voltage | Rated Current | Type      | Interrupt Rating  | Qty | Replacement period (year) |
|--|--------------|---------------|---------------|-----------|-------------------|-----|---------------------------|
| <b>Digital Acquisition Cabinet</b>         |              |               |               |           |                   |     |                           |
| Fuse, FLM30                                | 072-01659-88 | 250 V         | 30 A          | Slow Blow | 10,000 A@ 250 VAC | 2   | 2                         |
| Fuse, FLM10                                | 072-01659-83 | 250 V         | 10 A          | Slow Blow | 10,000 A@ 250 VAC | 2   | 2                         |
| Battery Unit, RBC6L                        | 074-81311-25 | -             | -             | -         | -                 | 1   | 2                         |
| Power Supply, SUA1500JB                    | 074-81311-09 | -             | -             | -         | -                 | 1   | 5                         |
| Battery Unit, BNB300S                      | 074-81330-71 | -             | -             | -         | -                 | 1   | 2                         |
| Power Supply, BN150S                       | 074-81330-03 | -             | -             | -         | -                 | 1   | 5                         |
| Battery, CR2032H                           | 074-73306-08 | -             | -             | -         | -                 | 2   | 2                         |
| HDD  | -            | -             | -             | -         | -                 | 1   | 5                         |
| IOF FAN2 ASSY F                            | 565-13125-01 | -             | -             | -         | -                 | 1   | 5                         |
| I-SYS FAN2 ASSY F                          | 565-13126-01 | -             | -             | -         | -                 | 1   | 5                         |
| <b>Digital Acquisition Cabinet Lateral</b> |              |               |               |           |                   |     |                           |
| Fuse, FLM30                                | 072-01659-88 | 250 V         | 30 A          | Slow Blow | 10,000 A@ 250 VAC | 2   | 2                         |
| Fuse, FLM10                                | 072-01659-83 | 250 V         | 10 A          | Slow Blow | 10,000 A@ 250 VAC | 2   | 2                         |
| IOF FAN2 ASSY L                            | 565-13125-02 | -             | -             | -         | -                 | 1   | 5                         |
| I-SYS FAN2 ASSY L                          | 565-13126-02 | -             | -             | -         | -                 | 1   | 5                         |
| <b>Digital Reference Cabinet</b>           |              |               |               |           |                   |     |                           |
| Fuse, 314 015                              | 072-01670-16 | 250 V         | 15 A          | Fast Blow | 750 A@ 250 VAC    | 1   | 2                         |
| Battery, CR2032H                           | 074-73306-08 | -             | -             | -         | -                 | 1   | 2                         |
| Mouse                                      | -            | -             | -             | -         | -                 | 1   | 5                         |

| Part                                       | Part No.     | Rated Voltage | Rated Current | Type      | Interrupt Rating  | Qty | Replacement period (year) |
|--|--------------|---------------|---------------|-----------|-------------------|-----|---------------------------|
| Keyboard                                   | -            | -             | -             | -         | -                 | 1   | 5                         |
| HDD  | -            | -             | -             | -         | -                 | 1   | 5                         |
| <b>Integrated Cabinet</b>                  |              |               |               |           |                   |     |                           |
| Fuse, FLM30                                | 072-01659-88 | 250 V         | 30 A          | Slow Blow | 10,000 A@ 250 VAC | 2   | 2                         |
| Fuse, FLM315                               | 072-01659-85 | 250 V         | 10 A          | Slow Blow | 10,000 A@ 250 VAC | 2   | 2                         |
| Fuse, FLM10                                | 072-01659-83 | 250 V         | 10 A          | Slow Blow | 10,000 A@ 250 VAC | 2   | 2                         |
| HDD  | -            | -             | -             | -         | -                 | 1   | 5                         |
| <b>IVR NEO</b>                             |              |               |               |           |                   |     |                           |
| WATERPROOF COVER NEO                       | 503-63823    | -             | -             | -         | -                 | 1   | 1                         |
| <b>FPD Cooling System</b>                  |              |               |               |           |                   |     |                           |
| Container for Anti-freezing solution, 10 L | 502-40190    | -             | -             | -         | -                 | 1   | 2                         |
| Container for Distilled Water, 1 L         | 502-40086    | -             | -             | -         | -                 | 1   | 2                         |
| <b>Transformer Box (Option)</b>            |              |               |               |           |                   |     |                           |
| Fuse, FLM30                                | 072-01659-88 | 250 V         | 10 A          | Slow Blow | 10,000A@ 250 VAC  | 1   | 2                         |
| Fuse, FLM20                                | 072-01659-86 | 250 V         | 10 A          | Slow Blow | 10,000A@ 250 VAC  | 4   | 2                         |
| Fuse, FLM10                                | 072-01659-83 | 250 V         | 10 A          | Slow Blow | 10,000 A@ 250 VAC | 2   | 2                         |
| <b>ECP</b>                                 |              |               |               |           |                   |     |                           |
| Fuse, 0326010.XP                           | 072-01665-34 | 250 V         | 10 A          | Slow Blow | 300 A@ 250 VAC    | 1   | 2                         |
| <b>Intercom (Option)</b>                   |              |               |               |           |                   |     |                           |
| Fuse, 313 002P                             | 072-01664-26 | 250 V         | 2 A           | Slow Blow | 100 A@ 250 VAC    | 2   | 2                         |



For IVR NEO and FPD cooling system, refer to the part number start with "56~" for RoHS.

This page is intentionally left blank.

# *Chapter 20*

---

## Appendix

### ■ Description

|                          |      |
|--------------------------|------|
| 20.1 DICOM Concepts..... | 20-2 |
| 20.2 Glossary.....       | 20-3 |

## 20.1 DICOM Concepts

This section provides basic DICOM conceptual information.

### 20.1.1 The DICOM Standard

DICOM is a communications and data encoding standard that provides interoperability between vendors. The DICOM standard specifies how to communicate data between stations rather than the type of data or how it is used.

To further understand DICOM, it is helpful to answer several key questions:

- **What is the DICOM standard?** The DICOM (Digital Imaging and Communications in Medicine) standard is a set of rules that allow medical images and associated information to be exchanged between imaging equipment, computers, and hospitals. The standard establishes a common language that permits medical images and information produced on one vendor's machine to be available for use on the digital system of another vendor.
- **Who wrote the DICOM Standard?** The DICOM Standard has evolved over the past ten years through meetings between medical imaging company representatives (represented through the National Electrical Manufacturers Association, NEMA) and representatives from major medical societies, including the American College of Cardiology, American College of Radiology, American Society of Echocardiography, European Society of Cardiology, and American Society of Cardiology.
- **Is compression of image data part of the standard?** Several compression schemes are defined in DICOM, but are limited by modality (X-ray, CAT, MRI, ultrasound). JPEG (Joint Photographic Experts Group) lossless compression, allows perfect reconstruction of the original image. It is mandated by the standard for cardiac angiography. JPEG lossy compression, which does not support perfect reconstruction, is a currently defined option only for echocardiography interchange; echocardiography also allows lossless data.
- **Can the exchange of digital images be standardized?** Concerned about incompatibility, members of the ACC, ASE, ESC and ASNC, working with the ACR-NEMA DICOM Standard committees began in 1992 to work on extensions to the standard to enable exchange of digital X-ray angiographic, echocardiography, and nuclear cardiac images. Each working group has recommended CD-R as an exchange medium because CD-R is non-erasable, sturdy, and easily obtained by both manufacturers and users. This allows each of these imaging modalities to be exchanged either separately or with each other, depending on the particular implementation of the standard.



Parts of this section are extracted from ACC's printed views on DICOM and the CD-R exchange media standards.

## 20.2 Glossary

All information in this glossary is provided in the context of medical imaging and this product.

|   |  |
|---|--|
| <b>ACC</b>                                | The American College of Cardiology.  |
| <b>ACR</b>                                | The American College of Radiology.   |
| <b>Cache</b>                              | A special reserved area on a hard disk that the product uses to temporarily store images and information.  |
| <b>Cath</b>                               | Cardiac Catheterization.   |
| <b>CAU</b>                                | Cranio-Caudal Projection.  |
| <b>CCD</b>                                | Charge-coupled device. A light-sensitive semiconductor device used in digital image capture.   |
| <b>CD</b><br><b>CD-ROM</b><br><b>CD-R</b> | Often used interchangeably, Compact Disc and Compact Disc Read only Memory can only be read and are produced in a pressing plant. Compact Disc Recordable (CD-R) discs can be written to in a computer drive but cannot be erased. |
| <b>Cine Run (loop)</b>                    | Video (multi-frame) images captured by medical imaging equipment.  |
| <b>CRA</b>                                | Caudo-Cranial Projection.  |
| <b>DA</b>                                 | Digital Angiography.   |
| <b>Diastole</b>                           | The dilation (expansion) of the heart chambers that occurs after contraction. Also referred to as “End Diastolic” state.   |
| <b>DICOM</b>                              | Digital Imaging and Communication in Medicine, is a set of rules that allow medical images and associated information to be exchanged between imaging equipment, computers, and hospitals.   |
| <b>DSA</b>                                | Digital Subtraction Angiography. A technique used to virtually remove constant structures of no diagnostic interest, enabling enhanced blood-vessel contrast.  |
| <b>ESC</b>                                | European Society of Cardiology.  |
| <b>FPD</b>                                | Flat Panel Detector.   |
| <b>Grayscale (Monochrome)</b>             | Black, white, and gray images that can have, for example, 256 unique shades with eight-bit data.   |
| <b>HIS</b>                                | Hospital information system.   |
| <b>IHE</b>                                | Integrating the Healthcare Enterprise.   |
| <b>JPEG</b>                               | Joint Photographic Experts Group, defines a set of standards for compression of still images and video (multi-frame) images.   |
| <b>LAO</b>                                | Left Anterior Oblique.   |
| <b>LIDM</b>                               | Live Image During Map  |
| <b>Lossless Compression</b>               | Digital compression in which data loss never occurs. Original data (image) condition is achieved after decompression.  |
| <b>Lossy Compression</b>                  | Digital compression which is not fully reversible but typically allows images to retain sufficient detail for analysis.  |

|                 |  |
|-----------------|--|
| <b>LUT</b>      | Look Up Table (used for brightness / contrast compensation).   |
| <b>LVA</b>      | Left Ventricular Analysis.   |
| <b>MR / MRN</b> | Medical Record / Medical Record Number.  |
| <b>MWM</b>      | Modality Worklist Management.  |
| <b>NEMA</b>     | National Electrical Manufacturers Association.   |
| <b>PPS</b>      | Performed Procedure Step.  |
| <b>QCA</b>      | Quantitative Coronary Analysis.  |
| <b>RAO</b>      | Right Anterior Oblique.  |
| <b>RDSR</b>     | DICOM Radiation Dose Structure Report  |
| <b>RSM-DSA</b>  | Real-time Smoothed Mask Digital Subtraction Angiography.   |
| <b>Stenosis</b> | An undesired narrowing of a patient's blood vessel.  |
| <b>Stent</b>    | A small tubular medical apparatus permanently inserted in a patient's blood vessel to prevent undesired narrowing.   |
| <b>Systole</b>  | The contraction of the heart chambers as blood is pumped into the aorta from the left ventricle and into the lungs from the right ventricle. Also referred to as "End Systolic" state. |