GPSL

GPSL-0497	R 1 x 100 mL + Std 1 x 5 mL
GPSL-0507	R 6 x 100 mL + Std 1 x 5 mL
GPSL-0707	R 4 x 250 mL + Std 1 x 5 mL
GPSL-0250	R 12 x 20 mL
GPSL-0455	R 6 x 45 mL
GPSL-0500	R 6 x 100 mL
GPSL-0700	R 4 x 250 mL

FTRO-GPSL-v20 (02/2022) (PIT-GPSL-4-v20)

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BIBLIOGRAFIE

1. Sacks, D.B., *Carbohydrates. Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A, Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 373.
2. Dods, R.F., *Diabetes Mellitus. Clinical Chemistry : Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds.), (2010), 729 and appendix
3. Wu, A. H. B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 444.
4. World Health Organization (WHO), *Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycemia*, (2006).
5. Trinder, P., *Ann. Clin. Biochem.*, (1969), 6, 24.
6. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
7. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), 59, 263.
8. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
9. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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0
PIT-GPSL

ADDENDUM

From the lot #22-0182 and beyond:
Use the version 23 of instructions for use
(page 2 to 5).

For the lot #22-0181 and previous:
keep using version 21 of instructions for use
(page 6 to 8).

PASL- 0230	R1 4 x 20 mL + R2 4 x 5 mL
PASL- 0400	R1 2 x 50 mL + R2 1 x 26 mL
PASL- 0420	R1 4 x 50 mL + R2 2 x 26 mL

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BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Dufour, D.R., *The Liver: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
3. *German Society for Clinical Chemistry, Z. Klin. Chem. Klin. Biochem.*, (1972), **10**, 281.
4. Wu, A.H.B., *Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders eds.), (2006), 78.
5. Rosalki, S.B., et al., *Clin. Chem.*, (1993), **39**, 648.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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FTRO-PASL-v21(12/2018)_PIT-PASL-4-v21

BIBLIOGRAFIE

- Henderson, A.R., Donald W.M., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 352.
- Scherwin, J.E., *Liver fonction. Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A, Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 492 and appendix.
- German Society for Clinical Chemistry, *Z. Klin. Chem. Klin. Biochem.*, (1972), **10**, 281.
- Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology, *Scand. J. Clin. Lab. Invest.*, (1974), **33**, 291.
- Guder, World Health Organization, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*, WHO/ DIL/LAB/99.1 Rev.2, (2002), 21.
- Rosalki, S.B., *Clin. Chem.*, (1993), **37**, 648.
- Vassault A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686.
- Vassault A., et al., *Ann. Biol. Clin.*, (1999), **57**, 685.
- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immuno-globulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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☞ BIBLIOGRAFIE

1. Rifai, N., Warnick, G.R., Remaley, A.T., *Lipids, lipoproteins, apolipoproteins and other cardiovascular risk factors*, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds.), (2008), 402.
2. Burnett, J.R., *Coronary Artery Disease: Lipid Metabolism. Clinical Chemistry: Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 691 and appendix.
3. Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III), *Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP)*, JAMA, (2001), **285**, 2486.
4. Fossati, P. & Prencipe, L., Clin. Chem., (1982), **28**, 2077.
5. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
6. Langlois, M.R., et al. for the European Atherosclerosis Society (EAS) and the European Federation of Clinical Chemistry and Laboratory Medicine (EFLM) Joint Consensus Initiative, Clin. Chem. Lab. Med., (2020), **58**, 496.
7. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
8. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press (1997).
9. Young D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press (1995).

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English - EN

INTENDED USE

This *in vitro* diagnostic device is intended for the calibration of ELITech Clinical Systems reagents, used for HDL or LDL cholesterol quantitative determination. This *in vitro* diagnostic device is for professional use only.

COMPOSITION

- Lyophilized product prepared from human serum.
- The HDL and LDL concentrations of this calibrator are lot-specific.

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- ELITech Clinical Systems reagent, used for HDL or LDL cholesterol quantitative determination (Ref. CHDL-XXXX or CLDL-XXXX).
- General Laboratory equipment.

TRACEABILITY

HDL LDL CALIBRATOR is traceable to the reference method approved by the CDC (Centers for Disease Control and Prevention).

PRECAUTIONS FOR USE AND WARNINGS

- Each unit of human blood used in the manufacture of these products was tested and found to be negative/non-reactive for the presence of HbsAg, HCV and HIV1/2. The methods used were FDA approved or CE compliant. Nevertheless, since the risk of infection cannot be fully excluded these products must be handled as potentially infectious. In case of exposure, follow the guidelines of the competent health authorities.
- Take precautions when handling broken glass vials as sharp edges can injure the user.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contaminations.
- Consult Safety Data Sheet (SDS) for a proper handling.

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PREPARATION

- Carefully open the vial, avoiding the loss of lyophilize.
- Add **exactly 1 mL of distilled or deionised water**.
- Carefully close the vial and dissolve the contents completely by occasional gentle stirring avoiding the formation of foam.
- Keep at room temperature for 30 minutes before use.

PRODUCT DETERIORATION

- The product may present a slightly hazy appearance after reconstitution. This has no effect on the performances of the product. All presence of particles would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter after reconstitution).
- Damage to the container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages).

STABILITY

Prior to reconstitution :

- Store at 2-8 °C and protect from light. Do not freeze.

- Do not use after the expiry date stated on the vial label.

After reconstitution :

- These products are stable 2 weeks at 2-8 °C or 10 weeks at -20°C.
- These products should be immediately and tightly capped to prevent contamination and evaporation.

PROCEDURE

To use HDL LDL CALIBRATOR, follow the procedure described in the instructions for use of the ELITech Clinical Systems reagent used.

LIMITATIONS

Using HDL LDL CALIBRATOR has been validated with the ELITech systems (Analyzers and reagents used for HDL or LDL cholesterol quantitative determination (Ref. CHDL-XXXX or CLDL-XXXX)).
Using any other system should be validated by the laboratory.

VALUES

The concentrations are indicated in the value sheet enclosed in the kit.
For users of Selectra instruments allowing data import for tests, calibrators and controls, use the corresponding barcode available on the value sheet.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the European Union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

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Holds Certificate No:

FM 703046

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

The design, manufacture, distribution and servicing of automated slide stainers, cytocentrifuges, cystic fibrosis sweat testing systems, and osmometers, and proprietary standards, controls disposables and reagents for use with these types of equipment. Manufacture and distribution of controls, standards, consumables, accessories and supplies for in vitro diagnostic systems, laboratory equipment, and erythrocyte sedimentation rate test systems.

For and on behalf of BSI:

Gary E Slack, Senior Vice President - Medical Devices

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BIBLIOGRAFIE

1. Wu, A.H.B., Tietz Clinical guide to laboratory test, 4th Ed., (W.B. Saunders Company), (2006), 648.
2. Sanhai, W.R., Eloff, B.C., Christenson, R.H., Cardiac and Muscle Disease, Clinical Chemistry: Theory Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby, Inc.), (2010), 677, appendix.
3. Panteghini, M., Bais, R., Enzyme, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (Saunders), (2008), 317.
4. Schumann, G. *et al.* Clin Chem Lab Med., (2002), **40**, 643.
5. Guder, W.G., *et al.*, Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. (2002). WHO/DIL/LAB/99.1 Rev.2
6. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), **23** (16).
7. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), **24** (34).
8. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
9. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
10. Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
11. Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), **59**, 263.
12. Young, D. S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).
13. Young, D. S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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Instrument Training

Vital Scientific BV hereby declares that the participant has attended a four days seminar for service engineers and the participant is now a certified engineer for the declared instruments.


Participant: Mr. A. Legun

Company: Global Biomarketing Group-Moldova SRL
Moldova

Instrument: Vitalab: XL Series
E Series
Junior Series
Dry ISE
Micro Series
ProXS


Date of training: April 20th – April 23rd, 2010

System Support Manager:

Jan Oostendorp

System Support Engineer:

Frank v.d. Kerput

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

Phosphorus is present throughout the human body as organic or inorganic phosphate: 80-85% in bones (hydroxyapatite), 15% in soft tissues (mainly in cells, integrated in nucleic acids, phospholipids, phospho-proteins, ATP, phosphorylated carbohydrates), less than 0.1% in extracellular fluids. In blood, phosphorus is mainly present as inorganic phosphate (Pi). Levels of Pi are closely linked to that of Calcium. An elevation of serum Pi can occur in renal failure, hypoparathyroidism or pseudo-hypoparathyroidism, excessive phosphate intake, leakage of intra-cellular phosphorus (e.g. untreated diabetic ketoacidosis) and vitamin D intoxication. A decrease of serum Pi can occur in hyperparathyroidism, renal tubule defects (Fanconi's syndrome), vitamin D deficiency (rickets, osteomalacia) and diabetic ketoacidosis induced by insulin.


Measuring Pi in urine can be recommended as an additional test to serum Pi. In practice, measurement of Pi in blood or urine is indicated to help diagnose pathologies involving abnormal serum calcium, for the monitoring of kidney pathologies, for the follow-up of diabetic patients or to assess gastrointestinal disorders.

LIMITATION OF USE

The quantitative assay of inorganic phosphorus alone can not be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

METHOD & PRINCIPLE ⁽⁴⁾

Phosphomolybdate - End point

COMPOSITION

Reagent: R	
Sulfuric acid	210 mmol/L
Ammonium molybdate	650 µmol/L
Surfactant	
☛ Standard: Std	
Phosphorus	4.7 mg/dL 1.52 mmol/L

☛MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- Analyzers or semi-automatic analyzers.
- General Laboratory equipment (e.g. pipette).
- Do not use materials that are not required as indicated above.

PRECAUTIONS FOR USE AND WARNINGS

- The reagent R is classified as hazardous.

- Consult Safety Data Sheet (SDS) for a proper handling.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.

STABILITY

Store at 2-25 °C and protect from light. Do not freeze.

Do not use after expiration dates indicated on the vial labels.

The standard should be immediately and tightly capped to prevent contamination and evaporation.

On board stability :

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagent and standard are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

SAMPLES

- Specimen ^(1,4)**
- Serum
 - Plasma (lithium heparin)
 - Urine
 - Using any other specimen type should be validated by the laboratory.
- Warnings and precautions**
- Collection of blood samples must be performed on fasting patients. ⁽³⁾
 - Separate blood samples from erythrocytes within 1 hour after collection. ⁽³⁾
 - Blood samples must be free from hemolysis. ⁽¹⁾
 - Urines must be acidified with HCl after collection (pH<3). ⁽³⁾
 - Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.
- Storage and stability ⁽²⁾**
- Serum/ plasma**
- 4 days at room temperature
 - 7 days at 2-8°C
 - Several months at -20°C
- Urine**
- 6 months

☛REFERENCE VALUES ⁽¹⁾

<i>Serum/plasma</i>	mg/dL	mmol/L
Adults	2.5 - 4.5	0.81 - 1.45
Children	4.0 - 7.0	1.29 - 2.26
<i>Urine (24h collection)</i>	g/24h	mmol/24h
	0.4 - 1.3	12.9 - 42.0
	For a urinary volume of 1.5 L/24h	
	mg/dL	mmol/L
	26.7 - 86.7	8.6 - 28.0

Note : *The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

☛PROCEDURE

Manual Procedure

Wavelength : 340 nm
Optical path : 1 cm
Sample/ Reagent ratio : 1:100
Temperature : 37 °C
Urine samples must be diluted 1:10 with NaCl 9 g/L solution before measurement
Read against reagent blank.

	CALIBRATION	TEST
Reagent R	1 000 µL	1 000 µL
Standard/ Calibrator	10 µL	-
Sample	-	10 µL

Mix and read the absorbances (A) after an incubation of 5 minutes.

Automatic Procedure

These reagents may be used on several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

Urine samples must be diluted 1:10 with NaCl 9 g/L solution before measurement. For users of Selectra TouchPro software, urine dilution is performed automatically.

☛INSTALLATION AND USE

Consult Selectra Pro operator manual.
Special Programming instructions: Programming special instructions is mandatory when some combinations of tests are performed together on the analyzer. Refer to instructions For Use of ACID SOLUTION & SYSTEM CLEANING SOLUTION for adequate programming (See PIT-SOL).

☛CALCULATION

A Sample x n n = Standard/Calibrator concentration
A Standard/ Calibrator

For the calculation of phosphorus concentration in urine, multiply the result by the dilution factor (10). For users of Selectra TouchPro software, the results take the dilution factor into account.

Conversion factor : mg/dL x 0.3229 = mmol/L

☛CALIBRATION

For the reference PHOS-0600 : ELICAL 2 and Phosphorus standard are traceable to SRM 956d reference material.
For the reference PHOS-0230 : ELICAL 2 is traceable to SRM 956d reference material.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay. Controls have to be performed :
- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements.
Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements (please refer to the Safety Data Sheet (SDS)).

PERFORMANCES

Performances were obtained on Selectra ProM, following CLSI technical recommendations, under controlled environmental conditions.

- Measuring range

a) Serum/Plasma

2.00 - 20.00 mg/dL (0.65 - 6.46 mmol/L).
Do not report results outside the measuring range.

b) Urine

10.0 - 200.0 mg/dL (3.2 - 64.6 mmol/L).
Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 1 000.0 mg/dL (322.9 mmol/L).

Do not report results outside this extended range.

☛For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

Limit of Detection (LoD) and Limit of Quantification (LoQ)

a) Serum/Plasma

LoD = 0.12 mg/dL (0.04 mmol/L)
LoQ = 1.00 mg/dL (0.32 mmol/L)

b) Urine

LoD = 0.2 mg/dL (0.1 mmol/L)
LoQ = 10.0 mg/dL (3.2 mmol/L)

- Precision

Imprecision data has been obtained on 2 Selectra ProM analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below :

a) Serum/Plasma

		Mean		Within-run	Total
	n	mg/dL	mmol/L	CV (%)	
Level 1	80	2.39	0.77	0.6	2.1
Level 2	80	4.77	1.54	0.4	2.3
Level 3	80	9.02	2.91	0.7	2.3

b) Urine

		Mean		Within-run	Total
	n	mg/dL	mmol/L	CV (%)	
Level 1	80	25.2	8.1	1.0	2.0
Level 2	80	72.5	23.4	0.6	1.0
Level 3	80	144.2	46.6	0.5	1.2

- Correlation

a) Serum/Plasma

A comparative study has been performed between PHOSPHORUS reagent on a Selectra ProM analyzer and a similar commercially available system on 100 human serum samples. The sample concentrations ranged from 2.01 to 18.82 mg/dL (0.65 - 6.08 mmol/L). The results are as follows :
Correlation coefficient : (r) = 0.999
Linear regression: y = 1.004 x + 0.05 mg/dL (0.02 mmol/L)

b) Urine

A comparative study has been performed between PHOSPHORUS reagent on a Selectra ProM analyzer and a similar commercially available system on 52 human urine samples. The sample concentrations ranged from 10.2 to 182.3 mg/dL (3.3 - 58.9 mmol/L). The results are as follows :
Correlation coefficient : (r) = 0.999
Linear regression: y = 0.883 x - 0.3 mg/dL (0.1 mmol/L)

- Limitations/Interferences

☛a) Serum/Plasma

- Hemolyzed samples should not be used since significant hemolysis may falsely increase the result due to the high level of phosphates in erythrocytes. ⁽¹⁾

- Studies have been performed to determine the level of interference from different compounds. The following inorganic phosphorus levels were tested: 2.17 mg/dL and 4.65 mg/dL.
- No significant interference is defined by a recovery ≤±10% of the initial value.
- **Unconjugated bilirubin**: No significant interference up to 15.0 mg/dL (257 µmol/L).
- **Conjugated bilirubin**: No significant interference up to 3.0 mg/dL (51 µmol/L)
- **Glucose** - No significant interference up to 528 mg/dL (29.3 mmol/L).
- **Triglycerides**: No significant interference up to 729 mg/dL (8.2 mmol/L).
- **Hemoglobin**: No significant interference up to 50 mg/dL.
- **Ascorbic acid**: No significant interference up to 20.0 mg/dL.
- **Methyl-Dopa**: No significant interference up to 1.0 mg/dL.

Do not use visibly turbid, icteric or hemolyzed samples.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results. ⁽⁵⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. ⁽⁶⁻⁷⁾

☛b) Urine

- Studies have been performed to determine the level of interference from different compounds. The following inorganic phosphorus levels were tested: 25.0 mg/dL and 150.0 mg/dL.
- No significant interference is defined by a recovery ≤±10% of the initial value.
- **Conjugated bilirubin**: No significant interference up to 22.1 mg/dL (378 µmol/L).
- **Hemoglobin**: No significant interference up to 500 mg/dL.
- **Uric acid** : No significant interference up to 120 mg/dL (7.1 mmol/L).
- **Urea** : No significant interference up to 6 000 mg/dL (999 mmol/L).
- **Ascorbic acid**: No significant interference up to 20.0 mg/dL.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. ⁽⁶⁻⁷⁾

- On board stability/Calibration frequency

On Board Stability: 28 days

Calibration frequency: 28 days

Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

These performances have been obtained using ELITech Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.

DECLARATION OF SERIOUS INCIDENT

Please notify the manufacturer (through your distributor) and competent authority of the Member State of the european union in which the user and/or the patient is established, of any serious incident that has occurred in relation to the device. For other jurisdictions, the declaration of serious incident should be in accordance with local, state and federal regulatory requirements. By reporting a serious incident, you provide information that can contribute to the safety of *in vitro* medical devices.

TECHNICAL ASSISTANCE

Contact your local distributor or ELITech Clinical Systems SAS (ccsupport@elitechgroup.com).

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BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFÍA/BIBLIOGRAFIA

1. Panteghini, M., Bais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders eds.), (2006), 64.
3. Dufour, R., *The liver: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc. eds.), (2010), 586 and appendix.
4. Schumann, G., *et al., Clin Chem Lab Med.*, (2002), **40**, 718.
5. Guder, W.G., *et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples.* (2002). WHO/DIL/ LAB/99.1 Rev.2.
6. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
7. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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English - EN

INTENDED USE
ELITech Clinical Systems AST/GOT 4+1 SL is an *in vitro* diagnostic reagent intended for the quantitative determination of aspartate aminotransferase (AST) in human serum and plasma samples on analyzers or semi-automatic analyzers.
This *in vitro* diagnostic device is for professional use only.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

Aspartate aminotransferase (AST) is a transaminase also known as glutamate oxaloacetate transaminase (GOT). AST catalyzes the transfer of the aminogroup of L-aspartate to α-ketoglutarate to give L-glutamate. AST is present in both mitochondria and cytoplasm and is widely distributed in the body, with high concentrations in the heart, liver, skeletal muscles and kidneys. Damage to any of these tissues may cause raised levels. In case of acute forms of hepatitis, especially viral hepatitis, the enzyme level is markedly elevated. In case of myocardial infarction, AST activity increases and reaches a peak after 18-24 hours. The activity falls back to normal after 4-5 days, provided no new infarct has occurred.



The following pathological states are examples of disorders also resulting in an increase of enzyme activity: other liver injury (e.g. chronic hepatitis, liver carcinoma, intake of alcohol, administration of drugs), muscular dystrophy, other damage to cardiac or skeletal muscles and acute pancreatitis.
AST measurement is indicated to screen patients at risk of liver diseases and to help diagnose or monitor liver diseases.

LIMITATION OF USE

The quantitative assay of aspartate aminotransferase (AST) alone can not be used to diagnose a disease or a specific pathology.
The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

METHOD & PRINCIPLE ⁽⁴⁾

IFCC method without pyridoxal phosphate (P-5'-P). Kinetic.

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PRECAUTIONS FOR USE AND WARNINGS

- Consult Safety Data Sheet (SDS) for a proper handling.
- The reagents contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.

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On board stability : The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

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SAMPLES

Specimen ^(2,5)

- Serum.
- Plasma (lithium heparin).
- Using any other specimen type should be validated by the laboratory.
- Warnings and precautions**
- Samples must be free from hemolysis. ^(1,3)
- Samples should be collected in accordance with Good Laboratory Practice and appropriate guidelines that may be in place.
- Storage and stability** ^(2,5)
- 24h at room temperature.
- 7 days at 2-8 °C.
- 3 month at -20 °C.

REFERENCE VALUES ⁽⁴⁾

<i>Serum/plasma</i>	U/L	µkat/L
Men	≤ 35	≤ 0.58
Women	≤ 31	≤ 0.52

Reference values for infants are higher than for adults.

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

Manual Procedure	
<i>Wavelelength :</i>	340 nm
<i>Optical path :</i>	1 cm
<i>Sample/reagent ratio :</i>	1:20
<i>Temperature:</i>	37 °C
<i>Read against distilled water.</i>	

Working reagent (4 volumes of R1 + 1 volume of R2)	1000 µL
Sample	50 µL

Mix and after 1 minute incubation, read absorbance at 1 minute intervals during 3 minutes. Calculate the change of absorbances per minute (ΔA/min)

Automatic Procedure

- These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

- High AST values may induce falsely low results due to the depletion of the substrate (total consumption of NADH before reading of the result). For ELITech Selectra Analyzers, the application contains a specific alarm to warn the users.

CALCULATION

Activity (U/L) = ΔA /min x 3333

Conversion factor : U/L x 0.0167 = µkat/L

CALIBRATION

ELICAL 2 is traceable to IFCC reference method.

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QUALITY CONTROL

It is recommended that quality control sera such as ELITROL I and ELITROL II be used to monitor the performance of the assay. Controls have to be performed :

- prior to assaying patient samples,
- at least once per day,
- after every calibration,
- and/or in accordance with laboratory and regulatory requirements.

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- Measuring range

10.0 to 450.0 U/L (0.17 to 7.50 µkat/L). Samples having greater concentrations should be diluted 1:10 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 4500.0 U/L (75.00 µkat/L).

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Representative results are presented in the following table.

		Mean		Within-run	Total
	n	U/L	µkat/L	CV (%)	
Level 1	80	34.1	0.57	1.7	3.4
Level 2	80	67.9	1.13	0.8	1.9
Level 3	80	353.6	5.89	0.4	2.0

- Correlation

A comparative study has been performed between AST/GOT 4+1 SL reagent on a Selectra ProM analyzer and a similar commercially available system on 114 human serum samples. The sample concentrations ranged from 10.0 to 529.8 U/L (0.17 - 8.83 µkat/L). The results are as follows :
Correlation coefficient : (r) = 0.999
Linear regression: y = 0.927 x - 0.3 U/L (0.01 µkat/L)

- Limitations/Interferences

- AST can be underestimated in case of severe vitamin B6 deficiency. ⁽³⁾

- Hemolyzed samples should not be used since significant hemolysis may increase AST concentration because of high levels of AST in erythrocytes. ⁽³⁾

- Studies have been performed to determine the level of interference from different compounds. The following AST levels were tested: 35.0 U/L and 350.0 U/L.

No significant interference is defined by a recovery ≤±10% of the initial value.
Unconjugated Bilirubin: No significant interference up to 30.0 mg/dL (513 µmol/L).
Conjugated Bilirubin: No significant interference up to 29.5 mg/dL (504 µmol/L).
Triglycerides: No significant interference up to 2400 mg/dL (27.1 mmol/L).

Pyruvate: No significant interference up to 3.0 mg/dL (341 µmol/L).

Ascorbic acid: No significant interference up to 20.0 mg/dL.

Acetylsalicylic acid: No significant interference up to 200 mg/dL.

Acetaminophen: No significant interference up to 30.0 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results. ⁽⁶⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. ⁽⁷⁻⁸⁾

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On Board Stability: 28 days
Calibration frequency: 28 days
Recalibrate when reagent lots change, when quality control results fall outside the established range and after a maintenance operation.

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BIBLIOGRAPHIE/BIBLIOGRAPHY

BIBLIOGRAFIA/BIBLIOGRAFIA

- 1.-Endres, D. B. & Rude, R.K., *Disorders of Bone*. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E. (W.B. Saunders eds), (2008), 711.
- 2.- Itani O. & Tsang, C., *Bone Disease. Clinical Chemistry: Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 614 and appendix.
- 3.-Swaminathan, R., *Clin. Biochem. Rev.*, (2003), 24, 47.
- 4.- Liamis, G. et al., Am. J. Nephrol., (2013), 38, 50.
- 5.- Mann, C.K, Yoe, J.H., *Spectrophotometric determination of magnesium with sodium 1-azo-2-hydroxy-3-(2,4-dimethylcarboxanilido)-naphthalene-1'-(2-hydroxybenzene-5-sulfonate)*, *Anal. Chem.*, (1956), 28, 202.
- 6.- Guder W.G., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
- 7.- Wu, A.H.B. *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds, Philadelphia USA), (2006), 706.
8. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), 59, 263.

9. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCC Press (1997).
10. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCC Press (1995)

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English - EN

INTENDED USE

ELITech Clinical Systems MAGNESIUM XB is an *in vitro* diagnostic reagent intended for the quantitative determination of magnesium in human serum, plasma and urine samples on analyzers or semi-automatic analyzers.

The standard is intended for the calibration of reagent. These *in vitro* diagnostic devices are for professional use only.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾

Magnesium is the fourth most abundant cation in the human body. It is the cofactor for many enzymatic systems including ATP-dependent enzymes and it also plays an active role in bone mineral homeostasis. Less than 1% of the total magnesium of the human body is carried by blood where it is mainly found as free ion, but also protein-bound (mainly albumin) or complexed with various anions. Magnesemia measures total magnesium but only free magnesium is biologically active. Hypomagnesemia can be due to important renal losses (excessive intake of diuretics), gastrointestinal disorders (malabsorption, diarrheas), endocrine disorders, less often due to inadequate intake. Hypomagnesemia is often associated with hypokalemia and/or hypocalcemia. Hypermagnesemia is rare and is caused almost always by excessive intake (parenteral therapy, excessive intake of magnesium-containing drugs) or renal insufficiency. The measure of urinary magnesium allows, in presence of magnesium deficiency, to identify a kidney etiology. In clinical practice, measurement of magnesium is indicated to help diagnose the causes of abnormal calcium and/or potassium levels or symptoms suggesting hypo or hypermagnesemia as well as for monitoring magnesium-and/or calcium-based treatments.

The measure of urinary magnesium allows, in presence of magnesium deficiency, to identify a kidney etiology. In clinical practice, measurement of magnesium is indicated to help diagnose the causes of abnormal calcium and/or potassium levels or symptoms suggesting hypo or hypermagnesemia as well as for monitoring magnesium-and/or calcium-based treatments.

LIMITATION OF USE

Because hypoalbuminemia may trigger pseudo-hypomagnesemia, serum magnesium results must be interpreted in regards of serum total protein and/or albumin levels. ⁽⁴⁾

The quantitative assay of MAGNESIUM XB alone can not be used to diagnose a disease or a specific pathology. The results must be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

METHOD & PRINCIPLE ⁽⁵⁾

Xylydyl Blue - End Point

Xylydyl blue in the reagent combines with the magnesium from the sample to form a red-purple chelate. EGTA is used to complex calcium and thus prevents it from interfering with the test.

The simultaneous increase in absorbance at 505-510 nm and decrease of the 620-630 nm absorbance are proportional to the magnesium concentration in the sample.

COMPOSITION

Reagent: R

AMP buffer, pH 11.2 120 µmol/L
Xylydyl Blue < 0.1 % (w/w)
Sodium azide
Also contains EGTA and surfactants for optimal performance

AMP: 2-Amino-2-methyl-1-propanol

Standard: Std
Magnesium 2.0 mg/dL
823 µmol/L

MATERIALS REQUIRED BUT NOT PROVIDED

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PRECAUTIONS FOR USE AND WARNINGS

- Reagent R contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.
- Consult Safety Data Sheet (SDS) for a proper handling.

STABILITY

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels.

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PRODUCT DETERIORATION

- These products should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of contamination or damage (e.g. particle matter).
- Damage to the product container may impact on product performance. Do not use the product if there is physical evidence of deterioration (e.g. leakages or punctured container).

Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take necessary corrective measures.

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a) *Serum/Plasma*

0.20 - 5.00 mg/dL (0.08 - 2.06 mmol/L)

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 25 mg/dL (10.28 mmol/L).

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For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

b) *Urine*

1.0 - 20.0 mg/dL (0.41-8.23 mmol/L)

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 100 mg/dL (41.14 mmol/L).

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For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

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LoQ = 1.0 mg/dL (0.41 mmol/L)

☛- Precision

Imprecision data has been obtained on 2 Selectra Pro analyzers over 20 days (2 runs per day, tests performed in duplicate).

Representative results are presented below.

a) *Serum/Plasma*

		Mean		Within-run	Total
	n	mg/dL	mmol/L	CV (%)	
Low level	80	1.55	0.64	1.0	3.5
Medium level	80	2.51	1.03	1.3	3.7
High level	80	3.84	1.58	0.7	3.0

b) *Urine*

		Mean		Within-run	Total
	n	mg/dL	mmol/L	CV (%)	
Low level	80	1.2	0.49	2.9	9.2
Medium level	80	5.0	2.06	0.7	4.0
High level	80	16.0	6.58	0.8	3.3

- Correlation

a) *Serum/Plasma*

A comparative study has been performed between MAGNESIUM XB reagent on a Selectra ProM analyzer and a similar commercially available system on 102 human serum samples.

The sample concentrations ranged from 0.27 to 4.99 mg/dL (0.11 - 2.05 mmol/L).

The results are as follows :

Correlation coefficient : (r) = 0.996

Linear regression: y = 1.041x - 0.01 mg/dL (0.00 mmol/L).

b) *Urine*

A comparative study has been performed between MAGNESIUM XB reagent on a Selectra ProM analyzer and a similar commercially available system on 80 human urine samples.

The sample concentrations ranged from 1.7 to 19.6 mg/dL (0.70 - 8.06 mmol/L).

The results are as follows :

Correlation coefficient : (r) = 0.999

Linear regression: y = 1.011x + 0.1 (0.00 mmol/L).

☛- Limitations/Interferences

a) *Serum/Plasma*

- Studies have been performed to determine the level of interference from different compounds. The following analyte levels were tested: 1.50 mg/dL and 3.90 mg/dL
No significant interference is defined by a recovery $\leq$ 10% of the initial value.
Triglycerides : No significant interference up to 3000 mg/dL(33.90 mmol/L).
Unconjugated bilirubin : No significant interference up to 30.0 mg/dL (513 µmol/L).
Conjugated bilirubin : No significant interference up to 29.5 mg/dL (504 µmol/L).
Hemoglobin : No significant interference up to 500 mg/dL.
Calcium : No significant interference up to 20.0 mg/dL (4.99 mmol/L).
Ascorbic acid : No significant interference up to 19.8 mg/dL.
Acetaminophen : No significant interference up to 200 mg/dL.

Acetylsalicylic acid : No significant interference up to 200 mg/dL

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results. ⁽⁶⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. ⁽⁹⁻¹⁰⁾

b) *Urine*

- Studies have been performed to determine the level of interference from different compounds. The following magnesium levels were tested: 1.2 mg/dL and 15 mg/dL.

No significant interference is defined by a recovery $\leq$ 10% of the initial value.
Conjugated bilirubin : Aucune interférence significative jusqu'à 29.5 mg/dL (504 µmol/L).
Hemoglobin : No significant interference up to 500 mg/dL.
Calcium : No significant interference up to 60.0 mg/dL (14.97 mmol/L).
Uric Acid : No significant interference up to100 mg/dL (5.95 mmol/L).
Urea : No significant interference up to 5000 mg/dL (832.50 mmol/L).
Ascorbic acid : No significant interference up to 19.8 mg/dL.
pH: No significant interference for pH values ranging between 2.5 à 6.0

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. ⁽⁹⁻¹⁰⁾

- On board stability/Calibration frequency
On Board Stability: 14 days
Calibration frequency: 7 days

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TECHNICAL ASSISTANCE

Contact your local distributor or ELITech Clinical Systems SAS (CCsupport@elitechgroup.com).

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To: Whom it May Concern

Regulatory status of parts & accessories

As mentioned on the current Declarations of Conformity of our Clinical Chemistry Analyzers also the accessories conform to the provisions of the EU Directive on In Vitro Diagnostic Medical Devices (98/79/EC). This applies to the parts and accessories as mentioned in the attached list.

'IVD accessory' means an article which, whilst not being an IVD medical device, is intended specifically by its manufacturer to be used together with an IVD device to enable that IVD device to be used in accordance with its intended purpose.

ELITechGroup B.V.

Adriaan P. Intveld
Manager Quality Assurance & Regulatory Affairs

Part number	Description	IVD medical device	IVD accessory	general laboratory use	spare part	supporting part
1540-001	Anti-Slip sheet					✓
2206-007	Cooling Liquid (1 L)					✓
3062-021	Sample cup (1000 pcs)		✓			
3062-033	Sample tube 6 ml (500 pcs)					✓
3062-040	Water container 10 L					✓
3062-041	Water container 5 L					✓
3066-155	Syringe 100 µl		✓			
3066-156	Syringe 1 ml		✓			
3069-040	Keyboard Dust cover					✓
3069-047	Keyboard Dust cover					✓
3070-518	Cap holder					✓
3070-538	Cap rotor Left					✓
3070-539	Cap rotor right					✓
3201-002	Dichromate 8 Abs (25ml)		✓			
3365-192	USB Stick					✓
3374-003	Mains cable (USA)					✓
3374-059	Pumpunit cable		✓			
3374-066	Mains cable					✓
3374-097	Serial Null-modem cable					✓
3374-286	USB Extension cable					✓
4804-038	Reagent identification Disc					✓
6001-826	Diluted Waste container		✓			
6001-827	Concentrated Waste container		✓			
6001-860	Water container		✓			
6001-861	Tube assy (analyser)		✓			
6001-872	Tube assy (cooling unit)		✓			
6002-102	Assorter unit				✓	
6002-386	System software on CD		✓			
6002-706	Reaction Rotor set (3 pcs)		✓			
6002-726	System Disc		✓			
6002-817	Bottle 30 ml (20 pcs)		✓			
6002-818	Bottle 15 ml (20 pcs)		✓			
6002-904	Water container 5 L		✓			
6002-910	Assorter unit				✓	
6002-913	External tubing		✓			
6003-074	System software on USB stick		✓			
6003-444	Diluted Waste Container 5 L		✓			
6003-466	Keyboard Support option					✓
6003-797	CW Waste Container 2 L		✓			
6003-808	Assorter unit				✓	

ISO 9001 -NF EN ISO 13485

R E A G E N T S

Zone Industrielle – 61500 SEES – France
Tél. : + 33 (0)2 33 81 21 00 / Fax : + 33 (0)2 33 28 77 51

TO WHOM TO BE CONCERNED

We, Seppim S.A.S., manufacturers of Elitech Clinical Systems reagents, having our factory at Zone Industrielle, 61500 Sées - France, confirm that our clinical reagents have been validated on Vital Scientific equipment. As such available Elitech Clinical Systems reagent applications for Vital Scientific instruments are CE-IVD compliant.

Reagents, other than Elitech Clinical Systems reagents, are not validated on Vital Scientific equipments, and we also can't know the impact of other reagents on Vital Scientific equipments.

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22 mai 2012

Signed on behalf of the manufacturer
Valérie GOURDON
Regulatory Affairs Manager
COMPANY SEPPIM S.A.S

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Environment : **Laboratory**

EN 61326-1:2006 : Equipment for measurement, control and laboratory use
EN 61326-2-6:2006 : Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 2-6: particular requirements – In vitro diagnostic (IVD) medical equipment, from which:

EN 55011:2007	: Emission - Class A
+A2:2007	
EN 61000-3-2:2006	: Limit for harmonic currents emissions
EN 61000-3-3:1995	: Limitation of voltage fluctuations and flicker
+A1:2001+A2:2005	
EN 61000-4-2:1995	: Electrostatic discharge (ESD) immunity
A1:1998+A2:2001	
EN 61000-4-3:2006	: Radiated Electro-Magnetic field immunity
+A1:2008	
EN 61000-4-4:2004	: Electrical fast transient (EFT) immunity
EN 61000-4-5:2006	: Surge transient immunity
EN 61000-4-6:2007	: Conducted Radio-Frequency disturbances immunity
EN 61000-4-8:1993	: Power frequency magnetic field immunity
+A1:2001	
EN 61000-4-11:2004	: Voltage dips and interruptions immunity

The undersigned declares that the described product meets the requirements of the mentioned standards, based on a non-recurrent examination. The test results lay down in our test reports with reference 2129388.0501-QUA/EMC and 2136226.0501-QUA/EMC.

KEMA Quality B.V.
 (Notified Body EMC)
 Arnhem, September 20, 2010

A.T. van der Meijden
 Certification Manager EMC

Certificate nr. **2136226.0551-QUA/EMC**

Integral publication of this certificate and adjoining reports is allowed.

KEMA Quality B.V. Utrechtseweg 310, 6812 AR Arnhem P.O. Box 5185, 6802 ED Arnhem The Netherlands
 T +31 26 3 56 20 00 F +31 26 3 52 58 00 www.kemaquality.com Registered Arnhem 09085396

GISL-0250	R1 8 x 20 mL + R2 8 x 5 mL
GISL-0400	R1 2 x 50 mL + R2 1 x 26 mL
GISL-0420	R1 4 x 50 mL + R2 2 x 26 mL

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☞BIBLIOGRAFIE

1. Panteghini, M., Bais, R., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A., Ashwood, E.R., Bruns, D.E., (W.B. Saunders eds.), (2008), 317.
2. Wu, A.H.B., *Tietz Clinical guide to laboratory test*, 4th Ed., (W.B. Saunders eds.), (2006), 470.
3. Dufour, R., *The Liver: Function and Chemical Pathology, Clinical Chemistry: Theory Analysis, Correlation*, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 586 and appendix.
4. Schumann, G., *et al.*, *Clin Chem Lab Med.*, (2002), **40**, 734
5. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations*, WHO/DIL/LAB/99.1 Rev.2, (2002).
6. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
- 7 Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
8. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995)

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☞ BIBLIOGRAFIE

1. Dufour D. R., *The Liver: Function and Chemical Pathology. Clinical Chemistry: Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A, Pesce, A.J., (Mosby Inc. eds), (2010), 586 and appendix.
2. Wu, A.H.B., *Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds.), (2006), 916.
3. Dumas, B.T., *et al.*, *Clin. Chem.*, (1981), **27**, 1642.
4. Guder, W.G., *et al.*, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*.(2002). WHO/DIL/LAB/99.1 Rev.2.
5. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
6. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed, AACC Press (1997).
7. Young D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed, AACC Press (1995).

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Contactați distribuitorul dumneavoastră local sau ELITech Clinical Systems SAS (ccsupport@elitechgroup.com).