

Clinical Summary

Single-centre randomized trial comparing conventional chemoembolization versus doxorubicin-loaded polyethylene glycol microspheres for early and intermediate stage hepatocellular carcinoma

Aleksandar Gjoreski et al, *European Journal of Cancer Prevention* 2020, DOI: 10.1097/CEJ.0000000000000623

Background

According to Barcelona Clinic Liver Cancer classification, transarterial chemoembolization (TACE) is preferred treatment for stage B and in certain cases for stage A hepatocellular carcinoma (HCC). Conventional TACE (c-TACE) and drug-eluting microspheres TACE (DEM-TACE) are available intra-arterial therapies.

Objective

The primary aim of this study was to compare the 12- and 24-month survival rates between the two arms. Secondary endpoints were comparison of intensity and duration of post-embolization syndrome (PES) after c-TACE and DEM-TACE and reporting of any serious adverse events after both methods, impact on liver functions, number of treatments and duration of hospital admission.

Methods

- A prospective, single-centre, randomized trial
- A total of 60 patients with unresectable HCC were included during a 36 months period and randomized one-to-one to undergo c-TACE or DEM-TACE
- Chemoembolization sessions were repeated “on demand” every 3-6 weeks until complete response.
- Follow-up of at least 24 months after the treatment or until death.
- 28 patients from the c-TACE group received a mixture of 50–100 mg of doxorubicin emulsified with 10–15 mL Lipiodol in a ratio 1:2.5, followed by administration of Contour 45–355 µm in diameter or HydroPearl 75–400 µm.
- 32 patients from the DEM-TACE group received LifePearl™ microspheres, 100–400 µm, volume of 2–4 mL preloaded with 50–100 mg of liquid doxorubicin. When needed, additional embolic agents for bland embolization were used, mostly polyethylene glycol particles, 75–200µm, (HydroPearl) at the discretion of the operator.

Results

Aleksandar Gjoreski et al study (n=60) reported:

- Overall Survival was different for the two arms but didn't reached statistical significance. groups, and showed:
 - DEM-TACE: at 12 months 90.2 and at 24 months 76.8%
 - cTACE: at 12 months 85.7 and at 24 months 63.6%
- One session was performed on 7.14% vs 15.63% of patients. This could represent a CR rate after one single intervention that is double for DEM-TACE vs cTACE
- CR/PR/SD rates were not reported

- **Safety:**

- o No significant difference in terms of adverse events was found.
- o Post Embolization Syndrome (PES) symptoms were slightly more severe after c-TACE, in particular elevated temperature (n. 23 -82.1% vs n. 13 40.6%) (P=0.001)
- o DEM-TACE required shorter in-hospital stay. The mean in-hospital stay was 2.37+1.4 vs 3.07+1.3 days for DEM-TACE vs cTACE. (statistically significant).

- **Limitations:**

- o The bigger limitation is that no hypothesis, neither sample size is described in the study methods, hence the non-significant difference in OS at 12 and 24 months could be due to an underestimate of the number of patients to enrol in order to show a real difference.
- o The authors didn't report the tumour response rates
- o Differences at baseline for nodules diameters: 5.82+2.4 for cTACE vs 6.23+2.6 for DEM-TACE
- o N. of pts with 1 intervention: 7.14% vs 15.63%; that translate in a CR rate after one single intervention that is double for DEM-TACE. This is not captured by the statistical analysis but could be clinically relevant if confirmed on further comparative studies.

CONCLUSION

Aleksandar Gjoreski et al concluded that this study did not demonstrate any statistically significant difference between c-TACE and DEM-TACE techniques in terms of 12-and 24-month survival rates. In real good candidates, both TACE methods are extremely effective and well tolerated with a great proportion of survival after 24 months. The only objective advantage of DEM-TACE over c-TACE is the shorter in-hospital stay after treatment.

Key Takeaways

- According to Gjoreski et al, DEM TACE-and cTACE procedures are effective and well tolerated treatments in selected HCC patients.
- Both procedures have comparable overall survival rates with a trend (NS) for higher rates for DEM-TACE.
- DEM-TACE allows a significant shorter in-hospital stay and a significant difference in one of the PES symptoms (fever)
- Link to the full publication: <https://pubmed.ncbi.nlm.nih.gov/33038087/>

Use/Indications

European Economic Area Indications for use

LifePearl™ microspheres are indicated for embolization of blood vessels supplying primary hypervascular tumours or metastases in the liver. Note: LifePearl™ microspheres can be loaded with chemotherapeutic drugs. When used for drug loading, drug loading should be done under a physician's direction, choice and responsibility, based on type and dose of drug most beneficial to the patient. LifePearl™ microspheres are compatible with doxorubicin, epirubicin, idarubicin and irinotecan. LifePearl™ microspheres can be drug loaded prior to embolization and then, as a secondary action, elute a local, controlled, and sustained dose to the targeted tumour sites after embolization. LifePearl™ microspheres are not available for sale in all countries. This information is provided only in respect to markets where this product is approved or cleared.

This literature summary is not a systemic review. It is only an example of LifePearl microspheres related literatures.

The use of the LifePearl™ devices in combination with drugs is not cleared or approved in the USA by the Food and Drug Administration. LifePearl™ microspheres are not approved in Canada.

Please consult the indication of use with the IFU supplemented with the product.

Please contact your Terumo local sales representative for more information. All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners. Refer to Instructions for Use for Contraindications, Warnings and Precautions.

HydroPearl™ microspheres are intended to occlude blood vessels for the Therapeutic and adjunctive Purposes. For complete list of indications, please refer to the instructions for use provided with the product. HydroPearl™ microspheres are not commercially available in all countries. Please contact your Terumo local sales representative for more information.

Please contact your Terumo local sales representative for more information. All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners. Refer to Instructions for Use for additional information. HydroPearl™ microspheres are approved for U.S. sales.

Lipiodol® (Ethiodized Oil) Injection is a registered trademark of Guerbet LLC.

CE0297. ©2021 MicroVention Europe

MicroVention Europe
30 bis, rue du Vieil Abreuvour
78100 Saint-Germain-en-Laye - France
Tel: +33 (0)1 39 21 77 46

Terumo Europe N.V.
Interleuvenlaan 40
3001 Leuven - Belgium
Tel: +32 16 38 12 11

 **TERUMO**
INTERVENTIONAL
SYSTEMS