CONTERM 19L 2000250  
CONTERM 36L 2000251  
CONTERM 52L 2000252  
CONTERM 80L 2000253  
CONTERM 150L 2000254

INCUBAT 19L 2000260  
INCUBAT 36L 2000261  
INCUBAT 52L 2000262  
INCUBAT 80L 2000263  
INCUBAT 150L 2000264

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## GENERAL INFORMATION


- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please advise the distributor immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) If you have any doubts or enquiries, please contact with your supplier or J.P. Selecta's technical service.
- 5) **IMPORTANT! J.P. SELECTA WILL NOT ACCEPT ANY MACHINE TO BE REPAIRED IF IT IS NOT DULY CLEANED.**
- 6) If any modification, elimination or lacking in maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.
- 7) Do not use the machine with liquids which can give off vapours capable of making explosive mixtures.

## PACKING LIST

The standard equipment consists of the following components:

	<i>Piece</i>	<i>Code</i>	<i>Qty.</i>
<b>Incubator 2000250/2000260</b>	Shelf runners	2000011	4
	Shelves	2000021	2
<b>Incubator 2000251/2000261</b>	Shelf runners	2000012	4
	Shelves	2000022	2
<b>Incubator 2000252/2000262</b>	Shelf runners	2000012	4
	Shelves	2000024	2
<b>Incubator 2000253/2000263</b>	Shelf runners	2000013	4
	Shelves	2000023	2
<b>Incubator 2000254/2000264</b>	Shelf runners	2000015	4
	Shelves	2000025	2
<b>Instructions manual</b>		80368	1

## ACCESSORIES

PROGRAMMABLE TIMER OF 24 HOURS 2000009  
 For continuous cycle with on-off switch levers that can be operated every 15 minutes.

## EQUIPMENT DESCRIPTION

**CONTERM** series are universal ovens for various applications, such as parts drying. The temperature is adjustable from room temperature + 5°C to 250°C.

**INCUBAT** incubators are designed for bacteriological culture. The heat distribution is produced by natural convection. The incubator is provided with a tempered glass inner door that allows the visualization without loss of temperature. Adjustable temperature from room temperature + 5°C to 80°C.

## TECHNICAL FEATURES

Voltage supply 115-230 V 50/60 Hz as indicated on the machine characteristics plate.

CONTERM MODELS		2000250	2000251	2000252	2000253	2000254
Power (W)		580	870	980	1150	1900
Capacity (L)		19	36	52	52	150
Temperature		40°C to 250°C				
Spatial variations*		±1°C to 150°C				
Stability*		±0.5°C to 150°C				
Thermometer error		±0.1°C				
Weight (Kg)		26	36	46	54	78
Inner dimensions (cm)	Height	30	40	33	50	50
	Width	25	30	47	40	60
	Depth	25	30	33	40	50
External dimensions (cm)	Height	50	60	53	70	70
	Width	60	65	82	74	95
	Depth	44	49	52	59	68

INCUBAT MODELS		2000260	2000261	2000262	2000263	2000264
Power (W)		170	225	250	300	525
Capacity (L)		19	36	52	52	150
Temperature		25°C to 80°C				
Spatial variations*		±0.5°C to 37°C				
Stability*		±0.3°C to 37°C				
Thermometer error		±0.1°C				
Weight (Kg)		26	36	46	54	78
Inner dimensions (cm)	Height	30	40	33	50	50
	Width	25	30	47	40	60
	Depth	25	30	33	40	50
External dimensions (cm)	Height	51	60	53	70	70
	Width	57	62	79	72	92
	Depth	49	54	57	64	74

\* Measured according to the measurement procedure of J.P. SELECTA, s.a.u.

## INSTALLATION

Put the incubator on a flat horizontal, level surface, leaving a space of about 10 cm at the back of the machine.

Do not install the incubator near sunny or heating places.

For an optimum function, the ambient temperature must not be over 30°C.

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Do not use the incubator to dry or for thermal treatment of products which can give off vapours capable of making an explosive mixture.

Be sure that the voltage supply is the same as the one indicated on the characteristics plate.

Do not use the equipment if it is not earthed.


If you change the plug bear in mind the following :

Blue cable: Neutral.  
 Brown cable: Phase.  
 Yellow/green cable: Earth.

## OPERATION

CONTROL PANEL:

1. General switch
2. Indicator screen
3. Stop button / Start cycle
4. Reduce Value key
5. Increase Value key
6. Adjustment (change between temperature / time)
7. Safety thermostat adjustment screw.

### STARTING UP

1. Press the start switch (1) and the screen (2) will light up at the same time.
2. Press keys (5)/(4) until you select the desired temperature value in the green display (the upper one).
3. Press key (6) if you want to select an operating time.
4. Press key (3) to start the cycle.
5. The heating symbol  will light and blink.
6. Press the key (3) again to stop the incubation process. Note: during operation, the temperature or time adjustment value can be changed.

### COMMENTS:

- \* Do not use the incubator at a lower temperature than ambient +5°C, because the thermostat might have difficulties in regulating correctly.
- \* When using the time function, the green display (the upper one) will indicate the remaining time. If you want to see the target temperature value, press key (6).
- \* Place the load in the incubator trying to leave spaces for the air circulation. Do not load the incubator more than 70% of its useful chamber volume.

### ADJUSTING THE SAFETY THERMOSTAT:

From the factory, the thermostat is adjusted to its maximum temperature.

To adjust the safety thermostat temperature, do as follows:

1. With the incubator connected to the mains, turn the safety thermostat screw (7) clockwise to the end of the scale, with the help of a flat screwdriver.
2. Select the desired safety temperature using the control keys. This temperature has to be minimum at 3°C above the working temperature.
3. With the door closed, turn on the incubator with key (3).
4. When the incubator is stabilized at the above selected temperature, turn the safety thermostat screw (7) slowly to the left until it clicks, and the red indicator lamp lights up.
5. Stop the incubator and leave it to cool with the door open during 15 minutes.

6. Select the working temperature by using the temperature selection keys (4) & (5) (minimum at 3°C below the desired safety temperature).
7. Close the door and turn the incubator on again with key (3). Once the temperature is stabilized, the incubator is ready to work.

**NOTE: Safety thermostat is especially useful when the user works habitually at the same temperature.**

## ERROR MESSAGES

P.FAIL Power failure. The stove continues with the process.

**Err** will appear in green on the top line of the display.

TSEG Safety thermostat is open.

SENS Open temperature probe. Control does not detect the temperature probe.

HEAT The heating element does not work. The equipment detects that the temperature does not increase.

## MAINTENANCE

Before removing the casing disconnect the equipment from the mains.

**The manipulation of the internal electronic circuits of the incubator by unauthorized personnel can cause irreparable damage. Take it to one of the J.P. SELECTA, s.a.u. authorized technical services.**

CLEANING:

For the cleaning of the different parts, use the following products:

Cleaning of stainless steel: Alcohol.

Cleaning of plastic: Alcohol with cotton duster.

## ADJUSTMENTS

The offset (A) and the °C/°F unit (B) can be adjusted.

### A. Adjust the probe reading correction (offset)

Press keys (4) and (5) continuously and at the same time (about 5 seconds). The **OFFs** message will appear.

Use keys (4) / (5) to set the actual temperature reading value. Press key (3) to confirm.

### B. Presentation unit: °C / °F

Press key (6). The message **Unit** will appear. Using keys (4) / (5), select the unit: °C / °F. Confirm with key (3).

## GUARANTEE

This product is guaranteed for one year. The guarantee does not cover damages caused by incorrect use or causes beyond the control of J.P. SELECTA, s.a.u.

Any manipulation of the equipment by personnel not authorized by J.P. SELECTA, s.a.u. automatically cancels the guarantee.

## SPARE PARTS

To guarantee the safety of the equipment, all spare parts must be purchased from J.P. SELECTA, s.a.u.

<i>Description</i>	<i>Code</i>				
<i>Model CONTERM</i>	2000250	2000251	2000252	2000253	2000254
<i>Schuko connector</i>	7171				
<i>Gasket</i>	21062				
<i>Heating element</i>	39122	39118	39120	39119	39121
<i>Thermostat</i>	38500				
<i>Glass door</i>	--	--	--	--	--
<i>Safety thermostat</i>			43030		

<i>Description</i>	<i>Code</i>				
<i>Model INCUBAT</i>	2000260	2000261	2000262	2000263	2000264
<i>Schuko connector</i>	7171				
<i>Gasket</i>	21254				
<i>Heating element</i>	39122	39118	39120	39119	39121
<i>Thermostat</i>	38500				
<i>Glass door</i>	47180	47181	47182	47183	47184
<i>Safety thermostat</i>			43153		