

MUTARS®



Total Knee MK
Surgical technique



implantcoast



Total Knee MK

Surgical technique

MUTARS® was developed in co-operation with
Prof. Dr. W. Winkelmann* (former director)
and Prof. Dr. G. Gosheger (director), Clinic and
Polyclinic for General Orthopedics and Tumororthopedics
at the University Hospital of Münster, Germany.
MUTARS® has been in successful clinical use since 1992.
**now Consultant for Orthopaedic Oncology,
Schönklinik Eilbek*

Table of contents

| | |
|----------------------------|----|
| Preoperative Planning..... | 5 |
| System Overview..... | 6 |
| Compatibility Matrix..... | 7 |
| Assembling Options..... | 8 |
| Surgical Technique..... | 9 |
| Implantats..... | 24 |
| Instruments..... | 28 |

Nota Bene: The described surgical technique is the suggested treatment for uncomplicated procedures. In the final analysis the preferred treatment is that which addresses the needs of the individual patient.

Copyright note: MUTARS®, GenuX®, ACS®, implavit® and implatan® are registered trademarks of implantcast GmbH. The use and copy of the content of this brochure are only allowed with prior permit given by the implantcast GmbH.

CAUTION: Federal law (USA) restricts this device to sale on or by the order of a physician.



MUTARS® Total Knee MK - the modular tumor system



characteristics

- full constrained knee system
- Mobile-Bearing and Fixed Bearing PE-inserts
- up to 20° rotation (Mobile-Bearing)
- cemented and cementless
- 2 femoral length
- 1 tibial sizes

MUTARS® Total Knee MK

Preoperative Planning

Preoperative planning and precise surgical techniques are mandatory for optimal results. The instructions and the procedure given in the surgical technique to the system must be adhered to. Familiarity with the recommended surgical technique and its careful application is essential to achieve the best possible outcome.

Prior to surgery a surgical planning with regard to the dimensions of the prosthetic model and the positioning of the implant components in the bone has to be carried out by the surgeon. For this purpose, x-ray templates are available from implantcast GmbH.




Figure A: MUTARS® Total Knee MK
Implant in A/P-View



Figure B: MUTARS® Total Knee MK
Implant in M/L-View



System Overview



MUTARS® femoral stem

Ø 11-17mm cemented
length: 120, 160, 200, 240mm
Ø 10-20mm cementless
length: 120 mm

MUTARS® Distal Femur M-O-M

length: 90, 100 mm

MUTARS® GenuX® MK

PE-Insert
size 2
MB and FB

MUTARS® MK Proximal Tibia

MUTARS® tibial stem

length: 120mm
Ø 11-15mm cemented




MUTARS® Total Knee MK

Compatibility Matrix

Distal Femur M-O-M



Distal Femur
90 mm
110 mm




PE-inserts MK
MB and FB



Size 2



Proximal Tibia MK





Assembling Options

(Length in mm)

| Components | | | | |
|-------------------------------|--------------|-----------------|-----------------|--------|
| femoral reconstruction length | distal femur | connecting part | extension piece | screw |
| 100 mm | 90 | - | - | 25 |
| 120 mm | 110 | - | - | 45 |
| 140 mm | 90 | - | 40 | 65 |
| 160 mm | 110 | - | 40 | 85 |
| 180 mm | 110 | - | 60 | 105 |
| 200 mm | 110 | - | 80 | 125 |
| 220 mm | 110 | 100 | - | 45+25 |
| 240 mm | 110 | - | 80+40 | 165 |
| 260 mm | 110 | 100 | 40 | 45+65 |
| 280 mm | 110 | 100 | 60 | 45+85 |
| 300 mm | 110 | 100 | 80 | 45+105 |
| 320 mm | 110 | 100 | 60+40 | 45+125 |

| Components | | | | |
|------------------------------|----------------|---------------------------------|-----------------|-------|
| tibial reconstruction length | prox. tibia MK | connecting part for prox. tibia | extension piece | screw |
| 115 mm | x | 105 | - | 25 |
| 135 mm | x | 125 | - | 45 |
| 155 mm | x | 105 | 40 | 65 |
| 175 mm | x | 105 | 60 | 85 |
| 195 mm | x | 105 | 80 | 105 |
| 215 mm | x | 125 | 80 | 125 |
| 235 mm | x | 125 | 40 + 60 | 145 |

Note: Please notice that the amount of implants and instruments send with an individual shipment may differ from the information in the catalogue information of this brochure. Please make sure, during the preoperatively planning, that all necessary implants and instruments are available for the surgery.

MUTARS® Total Knee MK

Distal femoral resection

Resect the tumour and measure the length of the explant.
The minimum bone resection should be 120mm (or 100 mm if the special Distal Femur 90 mm is used, available on demand). Remove the menisci.

Femoral preparation

Cemented fixation

Ream the femoral medullary cavity preferably up to a depth of 130 mm with a flexible reamer 2 mm larger than the preoperatively chosen femoral stem.

Cementless fixation

Ream the femoral medullary cavity preferably up to a depth of 130 mm with a flexible reamer that is 1,5 mm smaller than the preoperatively chosen femoral stem.

Chamfering of the medullary cavity

Prepare the femoral medullary cavity with the medullary cavity reamer cross-hole inside of the medullary cavity chamfer.

Note

If no flexible reamer in the clinic to be present, they are available on request!





Rasping of the femoral cavity

Cemented Implantation

Assemble the femur rasp of the appropriate size (see tables below), the sleeve and the slide hammer. Lock the rasp on the slide hammer by using the engineers' wrench.

Mark the anterior aspect of the femoral bone to meet the correct antecurvation of the femur

Note

If you want to prepare for a cemented stem with the femoral rasp, please use the rasp which is 2 mm larger than the preoperatively chosen cemented femoral stem.

That will provide a cement mantle of 1 mm thickness.

| Stem size | Rasp size |
|-----------|-----------|
| 11 mm | 13 mm |
| 13 mm | 15 mm |
| 15 mm | 17 mm |
| 17 mm | 19 mm |

Cementless Implantation

Use of cementless stems


Use the femur rasp, of the same size as the preoperatively chosen femoral stem.

| Stem size | Rasp size |
|-----------|-----------|
| 12 mm | 12 mm |
| 13 mm | 13 mm |
| 14 mm | 14 mm |
| 15 mm | 15 mm |
| 16 mm | 16 mm |
| 17 mm | 17 mm |
| 18 mm | 18 mm |
| 19 mm | 19 mm |
| 20 mm | 20 mm |

Note

It is recommended to clean the rasp of bone chips during the rasping.

Leave the femoral rasp in the bone for the trialing.



MUTARS® Total Knee MK

Tibial resection

Resect the tumour and measure the length of the explant.
The minimum bone resection should be 115mm.
Remove the menisci.

Tibial bone preparation

Use the medullary cavity reamer cross-hole to prepare the tibial bone.

Cemented fixation

Ream the tibial medullary cavity preferably up to a depth of 130 mm with a rigid reamer 2 mm larger than the size of the tibial stem.

Cementless fixation

Ream the femoral medullary cavity preferably up to a depth of 130 mm with a flexible reamer that is 1,5 mm smaller than the preoperatively chosen femoral stem.



Note

The use of a tibial rasp for a cemented stem is optional. Generally you can proceed with the trial reduction.



Cemented preparation

If you want to prepare for a cemented stem with the tibial rasp, please use the rasp which is 2 mm larger than the preoperatively chosen cemented tibial stem.

That will provide a cement mantle of 1 mm thickness. Use the 16mm rasp to prepare for the 15 mm stem.


| Stem size | Rasp size |
|-----------|-----------|
| 11mm | 13mm |
| 13mm | 15mm |
| 15mm | 16mm |

Cementless preparation

Choose the tibial rasp of the preoperatively planned size.

Assemble the tibial rasp of the appropriated size (see table below), the sleeve and the slide hammer. Lock the rasp on the slide hammer by using the engineers' wrench.

| Stem size | Rasp Size |
|-----------|-----------|
| 12mm | 12mm |
| 13mm | 13mm |
| 14mm | 14mm |
| 15mm | 15mm |
| 16mm | 16mm |



MUTARS® Total Knee MK

Although the tibial stem is not curved it is recommended to mark the anterior aspect of the tibial bone to assure that the rotation of the final stem is corresponding to the rotation of the rasp.

Rasp the medullary cavity with the chosen tibial rasp. A carefully use of the slap hammer is recommended.

To prevent fractures of the cortical bone it is helpful to fix a bone forceps around the tibial bone while rasping.

Note

It is recommended to clean the rasp of bone chips during the rasping.

Leave the tibial rasp in the bone for the trialing.





Trial reduction

Mount the trial Distal Femur and the possibly needed extension pieces (possible enlargement from 20 to 260 mm; see table page 8) to the top of the rasp.

Note

Please insert the cemented stem (without cement) or the trial stem for trialing purposes.

At that stage the use of a screw is optional, because the teeth mechanism gives the assembly a reasonable stability.



Insert the trial-Insert size 2 (MB or FB) in the trial Proximal Tibia MK.

MUTARS® Total Knee MK

Attach the trial connecting part for the Proximal Tibia (length: 105 mm or 125 mm) to the tibial rasp.

Note

Connect the trial Proximal Tibia MK and the connecting part. At that stage a bar screw is not needed, because the teeth mechanism provides reasonable stability.

Perform a trial reduction and check the joint stability and the rotational alignment.

Adjust the rotation if necessary. If the joint line could not be restored correctly, it might be necessary to change the length of the tibial reconstruction by a change of the connecting part, or adding of an extension piece in conjunction with an enlarged tibial bone resection.

If the correct position of the joint, the optimal reconstruction length as well as the stability and the rotational alignment is achieved, please remove all implant components.

Remove the tibial rasp, but consider the rotational alignment marks. Leave the Proximal Tibia and the PE-Insert combined.





MUTARS® Total Knee MK


Use the instrument for the locking mechanism to insert the GenuX® MK trial coupling into the femoral trial component.

The trial coupling and the instrument for the locking mechanism are assembled as shown **1**. Turn the coupling in a way that it falls into the sleeve of the instrument **2**.




Afterwards insert the coupling into the femoral box with the knee in flexion. For fixation turn the coupling via the instrument by 90° clockwise.

The coupling is positioned correct, when the peg falls out of the sleeve of the instrument. Remove the instrument afterwards.



Place the peg of the coupling in the opening of the tibial trial component in a way that the hole is directed towards anterior. Use the setting instrument for the locking mechanism or the assembling forceps as aid in positioning.

Insert the positioner from ventral into the hole of the tibial trial component to ensure the correct positioning of the peg.



Lock the coupling from ventral with the trial screw for coupling and the hex screwdriver 3.5mm.

MUTARS® Total Knee MK

Assamble the trial Proximal Tibia MK to the trial connecting part for proximal tibia by using the socket wrench and check the joint stability in flexion and extension.

Afterwards remove all trial components.



Implantation of the tibial / femoral stem

Impact the MUTARS® stem.

For a cemented implantation, insert the bone cement and use the cemented stem which is 2 mm smaller than the previously used reamer or rasp.

Remove all instruments, especially during the cement hardening to prevent bending moments.




Implantation of the femoral component

Combine the Distal Femur and possibly needed extension pieces with the femoral stem. Make sure that the correct rotation of the distal femur is achieved. Insert the bar screw of the correct length (see table on page 8).

Lock the bar screw with the MUTARS® swing wrench.
Secure the assembly with the engineers' wrench.

Lock the safety screw in the same way.




MUTARS® Total Knee MK

Use the instrument for the locking mechanism to insert the GenuX® MK coupling into the femoral component. The coupling and the instrument are assembled as shown **1**. Turn the coupling and the instrument in a way that the coupling falls into the sleeve of the instrument **2**.




Afterwards insert the coupling into the femoral box with the knee in flexion. For fixation of the coupling turn the coupling by 180° clockwise by use of the positioner as lever. The coupling is positioned correct, when the peg of the coupling falls out of the sleeve of the instrument. Remove the instrument for the locking mechanism afterwards.



Select the tibial insert of the desired type (Mobile-Bearing or Fixed-Bearing PE-insert) and place it to the tibial component. The size of the PE-insert corresponds to the size of the tibial component.


Place the peg of the coupling in the opening of the Proximale Tibia in a way that the hole is directed towards anterior. Use the setting instrument for the locking mechanism or the assembling forceps as aid in positioning.

Insert the positioner from ventral into the hole of the Proximal Tibia to ensure the correct positioning of the peg.





MUTARS® Total Knee MK



Combine the connecting part for proximal tibia with the femoral stem. Make sure that the correct rotation of the connecting part is achieved. Insert the bar screw of the correct length.

Lock the bar screw with the MUTARS® swing wrench.
Secure the assembly with the engineers' wrench.

Lock the safety screw in the same way.



Slide over the attachment tube. The trevira tube should be turned up inward on the end. If necessary cut the tube to the correct length




MUTARS® Total Knee MK

Explantation of the components

In case of an explantation please unlock the tibial counter screw and the locking screw from ventral first by use of the hex screwdriver 3.5mm.

Use the instrument for the locking mechanism to remove the coupling from the Distal Femur MOM component.




Subsequently the safety screw and the M10 bar screw are removed. The Distal Femur MOM can be easily removed from the stem now.

For the removal of the tibial components, please attach the tibial extractor to the slap hammer and insert it from above into the locking hole of the tibial component. Then lock the extractor with the locking tab of ventral. Subsequently the tibial component can be removed by the use of slap hammer.

Remove the two anterior screws from the proximal tibia by using the socket wrench and remove the proximal tibia.

Further loose the M10 screw by using the socket wrench and while counter with the engineers' wrench.

The connecting part for proximal tibia can be easily removed from the stem now.



Fixation on the attachment tube


Please fix the tube to the upper part of the Proximal Tibia by using the previously inserted 4 sutures.

Fold the tube to achieve a very close covering of the components. Insert additional sutures around the attachment pads of the implant components.



To reconstruct the extensor mechanism it is mandatory to perform a gastrocnemius muscle transfer. Release the muscle at its distal insertion. Suture the muscle to the anterior portion of the attachment tube.

Reinsert the extensor structures to the gastrocnemius muscle and the tube to restore a reasonable function of the joint.





Implants

***S:** For anti-infective treatment, silver coated implants are available.

***N:** For anti-allergic treatment, TiN coated implants are available.

***SN:** Implants with Silver and TiN coating!



MUTARS® GenuX® MK MB PE Insert

mat.: UHMW-PE acc. to ISO 5834-2

size

2 5721-0102



MUTARS® GenuX® MK FB PE Insert

mat.: UHMW-PE acc. to ISO 5834-2

size

2 5721-0202



MUTARS® MK proximal tibia *S *SN

mat.: implatan®; TiAl₆V₄ acc. to ISO 5832-3

5750-0005



MUTARS® connecting part for modular proximal Tibia *S

mat.: implatan®; TiAl₆V₄ acc. to ISO 5832-3

5750-0105 105 mm

5750-0125 125 mm



MUTARS® extension piece *S

mat.: implatan®; TiAl₆V₄ acc. to ISO 5832-3

5772-2504 40 mm

5772-2506 60 mm

5772-2508 80 mm

5772-2510 100 mm

MUTARS® Total Knee MK



MUTARS® connecting part *S

mat.: *implatan®*; $TiAl_6V_4$ acc. to DIN ISO 5832-3

5730-0100 100mm



MUTARS® GenuX® MK coupling *N

mat.: *implavit®*, *CoCrMo* acc. to ISO 5832-12

5720-1210




MUTARS® tibial stem cemented *N

mat.: *implavit®*; *CoCrMo* acc. to ISO 5832-4

5750-0511 11 mm

5750-0513 13 mm

5750-0515 15 mm



MUTARS® tibial stem cementless

mat.: *implatan®*; $TiAl_6V_4$ acc. to ISO 5832-3 with HA-coating acc. to ISO 13779-2

5750-1512 12 mm

5750-1513 13 mm

5750-1514 14 mm

5750-1515 15 mm

5750-1516 16 mm



MUTARS® Total Knee MK



MUTARS® tibial stem cemented with TiN-coating and HA-collar

mat.: *implavit®*; CoCrMo acc. to ISO 5832-4
with TiN and HA coating acc. to ISO 13779-2

| | |
|-----------|-------|
| 5759-1211 | 11 mm |
| 5759-1213 | 13 mm |
| 5759-1215 | 15 mm |



MUTARS® Distal Femur M-O-M, incl. safety screw *S *N *SN

mat.: *implavit®*; CoCrMo acc. to ISO 5832-4

| | | |
|-----------|--------|-------|
| 5720-0045 | 110 mm | left |
| 5720-0040 | 110 mm | right |
| 5720-0047 | 90 mm | left |
| 5720-0042 | 90 mm | right |



MUTARS® screw

mat.: *implatan®*; $TiAl_6V_4$ acc. to ISO 5832-3

| | |
|-----------|------------|
| 5792-1002 | M10x 25 mm |
| 5792-1004 | M10x 45 mm |
| 5792-1006 | M10x 65 mm |
| 5792-1008 | M10x 85 mm |
| 5792-1010 | M10x105 mm |
| 5792-1012 | M10x125 mm |
| 5792-1014 | M10x145 mm |
| 5792-1016 | M10x165 mm |
| 5792-1018 | M10x185 mm |
| 5792-1020 | M10x205 mm |
| 5792-1022 | M10x225 mm |



MUTARS® femoral stem cemented *N

mat.: *implavit®*; CoCrMo acc. to ISO 5832-4

| | |
|-----------|-------|
| 5760-0011 | 11 mm |
| 5760-0013 | 13 mm |
| 5760-0015 | 15 mm |
| 5760-0017 | 17 mm |

MUTARS® Total Knee MK



MUTARS® femoral stem cementless

mat.: implatan®, TiAl₆V₄ acc. to ISO 5832-3 with HA-coating acc. to ISO 13779-2

| | |
|-----------|----------|
| 5760-0111 | 11x120mm |
| 5760-0012 | 12x120mm |
| 5760-0113 | 13x120mm |
| 5760-0014 | 14x120mm |
| 5760-0115 | 15x120mm |
| 5760-0016 | 16x120mm |
| 5760-0117 | 17x120mm |
| 5760-0018 | 18x120mm |
| 5760-0019 | 19x120mm |
| 5760-0020 | 20x120mm |



MUTARS® attachment tube


mat.: Polyethylenterephthalat (PET)

| | |
|-----------|----------|
| 5900-0300 | Ø: 35 mm |
| 5900-0310 | Ø: 55 mm |




MUTARS® Total Knee MK


Instruments




7999-5821 MUTARS® distal femur MK basic container




7999-5834 MUTARS® proximal tibia MK basic container



7999-5836 MUTARS® proximal tibia MK trial container




7999-5835 MUTARS® proximal tibia MK drill container



7999-7701 MUTARS® trial container



7999-5745 MUTARS® patella container



7999-5712 MUTARS® basic container

MUTARS® Total Knee MK

MUTARS® distal femur MK basic container 7999-5821

resection check long
4220-0318

external alignment host
4223-0004

pin inserter 3,2 mm
4223-0006

ic T-handle Zimmer-Jakobs
4223-0023

slap hammer short
4223-0031

adapter M8 for slap hammer
4223-0032

external alignment rod 6x400mm
4223-0035

pin extractor
7512-0800

MUTARS® MK trial distal femur right/left

7720-0140 110mm right

7720-0142 90mm right

7720-0145 110mm left

7720-0147 90mm left

MUTARS® positioner for locking mechanism
7610-0003

MUTARS® instrument for locking mechanism
7720-1201

MUTARS® assembling forceps
7720-1202

GenuX® MK trial coupling
7720-1210

GenuX® MK trial screw for coupling
7720-1213

setting instrument for locking mechanism
7751-1200

ic pin adapter
4223-0421

ic adapter
4223-0022

MUTARS® proximal tibia MK basic container 7999-5834

MUTARS® patella drill guide
7350-0000

MUTARS® patella - clamp
7352-0001

extractor universal
7512-2026

hexagon screw driver short 3,5 mm
0280-1007

resection check
4223-0009

MUTARS® impact and extract sleeve
7230-0000

MUTARS® socket wrench
7420-0000

MUTARS® swing wrench
7411-0000

MUTARS® engineers' wrench SW 24
7490-0000

MUTARS® slide hammer
7220-0001

MUTARS® universal impactor
7210-0000

handle for intramedullary plug
7512-4001

MUTARS® rasp for tibial stem

7750-0312 12 mm

7750-0313 13 mm

7750-0314 14 mm

7750-0315 15 mm

7750-0316 16 mm



MUTARS® Total Knee MK

MUTARS® proximal tibia MK trial container 7999-5836

MUTARS® MK trial proximal tibia
7750-0005



MUTARS® trial screw for trial connecting
part 100mm (2 pieces)
7730-0106



MUTARS® trial extension piece for mod. prox. tibia
7750-0105 105mm
7750-0125 125mm



MUTARS® trial extension piece
7772-2504 40mm
7772-2506 60mm
7772-2508 80mm



MUTARS® trial bar screw
7792-1002 M10x25mm
7792-1004 M10x45mm
7792-1006 M10x65mm
7792-1008 M10x85mm
7792-1010 M10x105mm
7792-1012 M10x125mm
7792-1014 M10x145mm



GenuX® MK MB trial PE-insert 2
7721-0102



GenuX® MK FB trial PE-insert 2
7721-0202



MUTARS® proximal tibia MK drill container 7999-5835

MUTARS® rigid reamer

| | |
|-----------|---------|
| 7700-2110 | 10,0 mm |
| 7700-2210 | 10,5 mm |
| 7700-2111 | 11,0 mm |
| 7700-2211 | 11,5 mm |
| 7700-2112 | 12,0 mm |
| 7700-2212 | 12,5 mm |
| 7700-2113 | 13,0 mm |
| 7700-2213 | 13,5 mm |
| 7700-2114 | 14,0 mm |
| 7700-2214 | 14,5 mm |
| 7700-2115 | 15,0 mm |
| 7700-2116 | 16,0 mm |
| 7700-2117 | 17,0 mm |



MUTARS® medullary cavity reamer
7760-0501



MUTARS® patella container 7999-5745

MUTARS® patella drill guide
7350-0000



MUTARS® patella - clamp
7352-0001



MUTARS® patella drill
7351-0000



MUTARS® Total Knee MK

MUTARS® trial container 7999-7701

MUTARS® trial prox. femur

7710-0205 50mm

7710-0207 70mm



MUTARS® trial reducer

7730-0220 20mm

7730-0230 30mm



MUTARS® trial connecting part 100 mm

7730-0100



MUTARS® trial extension piece for mod. prox. tibia

7750-0105 105mm

7750-0125 125mm



MUTARS® trial extension piece

7772-2504 40mm

7772-2506 60mm

7772-2508 80mm

7772-2510 100mm



MUTARS® trial femoral stem

7760-0011 11mm

7760-0013 13mm

7760-0015 15mm

7760-0017 17mm



MUTARS® trial bar screw

7792-1002 M10x25mm

7792-1004 M10x45mm

7792-1006 M10x65mm

7792-1008 M10x85mm

7792-1010 M10x105mm

7792-1012 M10x125mm

7792-1014 M10x145mm

7792-1016 M10x165mm

7792-1018 M10x185mm

7792-1020 M10x205mm



MUTARS® basic container 7999-5712

MUTARS® universal impactor

7210-0000



MUTARS® impact and extract sleeve

7230-0000



MUTARS® socket wrench

7420-0000



MUTARS® swing wrench

7411-0000



MUTARS® engineers' wrench SW 24

7490-0000



MUTARS® slide hammer

7220-0001



MUTARS® rasp for femoral stem

7760-0112 12 mm

7760-0113 13 mm

7760-0114 14 mm

7760-0115 15 mm

7760-0116 16 mm

7760-0117 17 mm

7760-0118 18 mm

7760-0119 19 mm

7760-0120 20 mm



handle for intramedullary plug

7512-4001



MUTARS® medullary cavity reamer cross-hole

4220-0000





implantcast GmbH
Lüneburger Schanze 26
21614 Buxtehude
phone.: +49 4161 744-0
fax: +49 4161 744-200
e-mail: info@implantcast.de
internet: www.implantcast.de



Your local distributor:

TKN/KOPE-200716

