DECLARATION DE CONFORMITE CE

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DECLARATION OF EC CONFORMITY

We, ELITech Clinical Systems SAS, Zone Industrielle 61500 SEES France, hereby certify, under our own responsibility, that the reagents belonging to Group 1 "MISCELLANEOUS METABOLITES", such as listed hereto, conform to the essential requirements of appendices I and III of European Directive 98/79/EC, relating to *in vitro* diagnostic medical devices and to the public health code.

This declaration is based on the contents of each DOS-CE-XXXX technical file and is supported by the certification of our quality system according to the standard NF EN ISO 13485 : 2016 (Certification valid until July 27th, 2020).

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BIBLIOGRAFIE

1. Scherwin, J.E., *Liver function. Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 492 and appendix.
2. Tietz, N.W., *Clinical guide to laboratory tests*, 3rd Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 518.
3. Doumas, B.T., et al., *Clin. Chem.*, (1981), **27**, 1642.
4. Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2.
5. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI (NCCLS) document EP6-A (2003), **23** (16).
6. *Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
7. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
8. *Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
9. *Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
10. Vassault, A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686.
11. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), **59**, 263.
12. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
13. Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

1. Sacks, D.B., *Carbohydrates*. Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2008), 373.
2. Dods, R.F., *Diabetes Mellitus*. Clinical Chemistry : Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A, Pesce, A.J., S.C., (Mosby Inc. eds St Louis USA), (2010), 729 and appendix
3. Wu, H.B., *General Clinical Tests*. Tietz Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 444.
4. Trinder, P., *Determination of glucose in blood using glucose oxidase with an -alternative oxygen -acceptor*. Ann. Clin. Biochem., (1969), 6, 24.
5. Guder W.G, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma, and serum samples*, World Health Organization , WHO/DIL/LAB/99.1 Rev.2, (2002).
6. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), 23 (16).
7. *Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), 24 (34).
8. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), 24 (25).
9. *Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), 22 (19).
10. *Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), 25 (27).
11. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), 59, 263.
12. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
13. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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BIBLIOGRAFIE

- Higgins, T., et al., *Hemoglobin, Iron, and Bilirubin*, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds, Philadelphia USA), (2008), 509.
- Dufour, D.R. *The liver: Function and chemical Pathology*, Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc.), (2010), 586 and appendix.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
- Wu, H.B., *General Clinical Tests*. Tietz Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 172.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).8.
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, Acta Clin Belg., (2004), **59**, 263.
- Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd edition, AACC Press (1997).
- Young D.S., Effects of drugs on clinical laboratory tests, 4th edition, AACC Press (1995).

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**GOST R CERTIFICATION SYSTEM
FEDERAL AGENCY FOR TECHNIQUE REGULATION AND
METROLOGY**

**VOLUNTARY CERTIFICATION SYSTEM
"SMK-STANDARD"**

Reg. No. POCC RU.31060.04ЖЖЮ0

Certification authority:

REG No. SMK STANDART.RU.0005

INTERNATIONAL CERTIFICATION CENTER Limited Liability Company

Address: 138, Naberezhnaya Obvodnogo Canala, block 1, office 421, St. Petersburg, 190020

phone: +7 (812) 438-76-71 standart@iso-smk.ru

Check the authenticity of the certificate in the register on the website <http://www.iso-smk.ru>

CERTIFICATE OF CONFORMITY

No. ST.RU.0001.M0013380

This Certificate of Conformity is issued to

Agat-Med, Ltd.

Address: 6, Glavnaya st., Moscow, 105173, Russia

TRN 7719187311 OGRN 1037739078970

Date of issue: 26.01.2018

Period of validity: 26.01.2021

This certificate certifies that:

Medical devices. Quality management system. System requirements for regulatory purposes in relation to the works in accordance with Annex 1 to this certificate

(the attachment is an integral part of the certificate)

CORRESPONDS TO THE REQUIREMENTS OF GOST ISO 13485-2011 (EN ISO 13485:2003)

Manager of "Expert" authority

V.V. Koptsev

O.V. Gundareva

THIS CERTIFICATE BINDS THE ORGANIZATION TO MAINTAIN THE WORKS PERFORMED ACCORDING TO THE STANDARD MENTIONED ABOVE, TO BE CONTROLLED BY THE INSTITUTIONAL BODY OF VOLUNTARY CERTIFICATION SYSTEM "SMK STANDARD" AND IS CONFIRMED DURING THE ANNUAL INSPECTION CONTROL

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Test strips for rapid determination of protein, ascorbic acid, glucose, pH-value in urine

Use

Screening test for detection of diabetes and diseases of kidney and urinary tract.
Only for use by qualified personnel.

Instructions for use

Dip the test strip for approximately 1 second into the fresh urine. Draw it across the rim of the container to remove excess urine. After 30 to 60 seconds compare the test strip with the color scale. The best time for comparison is after 30 seconds. Color changes that take place after more than 2 minutes are of no significance. When tested the urine should not be older than 2 hours.

Principle

Protein: The test is based on the "protein error" principle of indicators. The test zone is buffered to a constant pH value and changes color from yellow to greenish blue in the presence of albumin. Other proteins are indicated with less sensitivity.

Ascorbic acid: The detection is based on the decoloration of Tillmans reagent. In the presence of ascorbic acid a color change takes place from blue to red.

Glucose: The detection is based on the glucoseoxidase-peroxidase-chromogen reaction. Apart from glucose, no other compound in urine is known to give a positive reaction.

pH: The test paper contains indicators which clearly change color between pH 5 and pH 9 (from orange to green to turquoise).

Evaluation – Sources of Error

Protein: The minimum sensitivity of the test strip is 10 mg protein/dL urine. The color fields correspond to the following ranges of albumin concentrations:

negative, 30, 100, 500 mg/dL or
negative, 0,3, 1,0, 5,0 g/L

Falsely positive results are possible in alkaline urine samples (pH > 9), after infusions with polyvinylpyrrolidone (blood substitute), after intake of medicaments containing quinine and also by disinfectant residues in the urine sampling vessel. The protein coloration may be masked by the presence of medical dyes (e.g. methylene blue) or beetroot pigments.

Ascorbic acid: The color fields correspond to the following values:

0 (negative), 10(+), 20(++), 100 mg/dL or
0 (negative), 0,6(+), 1,1(++), 11 mmol/L

Only for information!

Glucose: Pathological glucose concentrations are indicated by a color change from green to bluish green. Yellow or greenish test fields should be considered negative or normal. The color fields correspond to the following ranges of glucose concentrations:

neg. (yellow), neg. or normal (greenish), 50, 150, 500, ≥ 1000 mg/dL or
neg. (yellow), neg. or normal (greenish), 2,8, 8,3, 27,8, ≥ 55,5 mmol/L

The influence of ascorbic acid (vitamin C) has been largely eliminated. An inhibitory effect is produced by gentisic acid. Falsely positive reactions can be produced by a residue of peroxide containing cleansing agents.

pH: The pH value of fresh urine of healthy people varies between pH 5 and pH 6. The color scale gives a clear distinction of pH value between pH 5 and pH 9.

Reacting substances

(Quantity resp. activity/cm² at time of impregnation)

Protein:	Glucose:	pH:	Methyl red	3 µg
tetrabromophenol blue	10 µg	7 U	methyl red	10 µg
Ascorbic acid:	1 U	96 µg		
2,6-dichlorophenolindophenol	7 µg			

Quality Control in professional use

The performance of the test strips should be confirmed by use of positive and negative control solutions. Positive and negative controls should be analyzed once a day, whenever a new bottle of strips is opened, whenever a new lot of strips is started, and every 30 days to check storage conditions. Each laboratory should establish its own goals for adequate standards of performance, and should question handling and testing procedures if these standards are not met.

Directions

In any case, in order to establish a final diagnosis and prescribe an appropriate therapy, the results obtained with test strips should be verified with other medical results.

The effect of medicaments or their metabolic products on the test is not known in all cases. In case of doubt it is recommended not to take the medicaments and then repeat the test.

Only use well washed and clean vessels for urine collection. The presence of usual urine preservatives will not affect the test results.

Remove only as many test strips as are required, and reseal the container immediately after use. Do not touch the test pads. Avoid exposing the strips to sunlight and moisture. Store the container below +30 °C in a dry place. The test strips are stable, when stored properly up to the date of expiry indicated.

The caps contain a non-poisonous and harmless desiccant. In case this desiccant is swallowed accidentally, then drink plenty of water.

Explanation of symbols can be found in the package insert.

Disposal: Please dispose all used test strips in accordance with your local laws and regulations.

Package units: Tubes of 50 and 100 test strips

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Tel. : 062 388 55 00 · Fax : 062 388 55 05 · sales-ch@mn-net.com

EC Certificate

Directive 98/79/EC Annex IV, excluding Sections 4 and 6
Full Quality Assurance System
In Vitro Diagnostic Medical Devices

Registration No.: HL 60119814 0001

Report No.: 21265422 001

Manufacturer: Macherey-Nagel GmbH & Co. KG
Neumann-Neander-Str. 6-8
52355 Düren
Deutschland

Products: Products for self-testing
(see attachment for products and sites included)
Replaces Certificate, Registration No.: HL 60076687 0001

Expiry Date: 2022-05-28

The Notified Body hereby declares that the requirements of Annex IV, excluding section 4 and 6 of the directive 98/79/EC have been met for the listed products. The above named manufacturer has established and applies a quality assurance system, which is subject to periodic surveillance, defined by Annex IV, section 5 of the aforementioned directive. For placing on the market of List A devices covered by this certificate an EC design-examination certificate according to Annex IV, section 4 and a verification of manufactured products according to section 6 is required.

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TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 98/79/EC concerning in vitro diagnostic medical devices with the identification number 0197.

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Products for self-testing:

- Single and multi-parameter disposable test strips for urine analysis
- indicator test strips and papers for measurement of pH in urine

Additional site for warehousing and logistics:

Bahnstr. 120
52355 Düren, Germany

Date: 2017-05-29

Notified Body

Dipl.-Ing. Sven Hoffmann

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TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg
Manufacturer: Macheray-Nagel GmbH & Co. KG