



Artis one

Maintenance Instructions – System
AX System Checkup and Maintenance

Standards

This document is IEC 62353:2014-09 compliant.

Document Version

Siemens Healthineers AG reserves the right to change its products and services at any time.

In addition, manuals are subject to change without notice. The hard copy documents correspond to the version at the time of system delivery and/or printout. Versions to hard copy documentation are not automatically distributed.

Please contact your local Siemens Healthineers office to order a current version or refer to our website ↗ <https://www.siemens-healthineers.com>.

Disclaimer

Siemens Healthineers AG provides this documentation “as is” without the assumption of any liability under any theory of law.

The content described herein shall be used by qualified personnel who are employed by Siemens Healthineers AG or one of its affiliates or who are otherwise authorized by Siemens Healthineers AG or its affiliates to use such documents.

Assemblers and other persons who are not employed by or otherwise directly affiliated with or authorized by Siemens Healthineers AG or one of its affiliates are not entitled to use this documentation without prior written authority.

This documentation may contain pictorial representations of work steps (pictures, videos, animations, etc.). Such pictorial representation is to be understood as complementing the documentation and is not to be used standalone. The pictorial representations are protected by copyright and it is not permitted to extract such content from the context of this documentation or use it on its own or distribute it independently.

Copyright

“© Siemens Healthineers, 2015” refers to the copyright of a Siemens Healthineers entity such as:

Siemens Healthineers AG - Germany

Siemens Shenzhen Magnetic Resonance Ltd. - China

Siemens Shanghai Medical Equipment Ltd. - China

Siemens Healthcare Private Ltd. - India

Siemens Medical Solutions USA Inc. - USA

Siemens Healthcare Diagnostics Inc. - USA

Siemens Healthcare Diagnostics Products GmbH - Germany

Corindus Inc. - USA

1	General Information	7
1.1	Prerequisites	7
1.2	General Information	7
1.2.1	General Safety Information	7
1.2.2	Manufacturer's Note	7
1.2.3	Maintenance Intervals	7
1.2.4	Notes regarding Maintenance	8
1.2.5	Explanation of Abbreviations in the Maintenance Instructions	8
1.3	Product-specific Information	8
1.3.1	Explanation of Abbreviations in the Maintenance Certificate	8
1.3.2	Maximum torque in Nm (table)	9
1.3.3	Tools, Aids, Lubricants, Cleaning Agents and Paints	9
1.3.4	Components Taken into Consideration	11
1.3.5	DHHS Regulations	12
1.3.6	Components affected by DHHS regulations	12
1.3.7	DHHS labels on the system	12
1.3.8	Corrective Measures/Adjustments	16
1.3.9	Cover Panels	17
1.4	Product-specific Safety Information	17
1.4.1	Patient Table	23
2	System Checkup	26
2.1	Special Requirements	26
2.1.1	Required Documents	26
2.1.2	Required Tools, Measuring Equipment and Aids for All Components	26
2.1.3	Required Materials for all Components	27
2.2	General Visual Check	27
2.2.1	Visual Check per IEC 62353	27
2.3	Image System (AXIS Smart)	28
2.3.1	Product-specific Safety Information	28
2.3.2	Mechanical Safety	29
2.4	Panoramic Display (option)	29
2.4.1	Product-specific Safety Information	29
2.4.2	Mechanical Safety	29
2.4.3	Electrical Safety	29
2.5	IOB/CCOM (RTC), Input/Output Box, Central Computer	30
2.5.1	Product-specific Safety Information	30
2.5.2	Electrical Safety	30
2.6	TFT Color Monitor (Examination Room)	30
2.6.1	Product-specific Safety Information	30
2.6.2	Mechanical Safety	31
2.7	21" Display (Control Room)	31
2.7.1	Mechanical Safety	31
2.8	X-ray Generator	31
2.8.1	Product-specific Safety Information	31
2.8.2	Electrical Safety	32
2.9	Tube Cooling Unit	33
2.9.1	Product-specific Safety Information	33
2.9.2	Visual check	34
2.10	DCS	34
2.10.1	Exterior Visual Inspection	34
2.10.2	Mechanical Safety	35

2.11	Upper Body Radiation Shield, Option	35
2.11.1	Product-specific Safety Information	35
2.11.2	Mechanical Safety	35
2.12	Stand	36
2.12.1	Product-specific Safety Information	36
2.12.2	Mechanical Safety	36
2.12.3	Cleaning and Lubricating the Rotation Drive	40
2.12.4	Cleaning and Lubricating the FD Receptor	41
2.12.5	Electrical Safety, Safety Devices	42
2.13	Patient Table	45
2.13.1	Product-specific Safety Information	45
2.13.2	Mechanical Safety	45
2.13.3	Electrical Safety	46
2.13.4	Preventive Lubrication	47
2.14	Final Work Steps	50
2.14.1	Installing the Cover Panels	50
2.14.2	Electrical Safety	51
3	Full Maintenance	54
3.1	Part 1 AXIS Smart/User Interfaces/IOB/Displays	55
3.1.1	Special Requirements	55
3.1.2	General Visual Check	56
3.1.3	Image System (AXIS Smart)	57
3.1.4	Panoramic Display (option)	60
3.1.5	User Interfaces	62
3.1.6	IOB/CCOM (RTC), Input/Output Box, Central Computer	67
3.1.7	MCI	68
3.1.8	TFT Color Monitor (Examination Room)	68
3.1.9	21" Display (Control Room)	71
3.1.10	Final Work Steps	72
3.2	Part 2: Generator/DCS/Upper Body Radiation Shield	72
3.2.1	Special Requirements	72
3.2.2	X-ray generator	74
3.2.3	Tube Cooling Unit	83
3.2.4	DCS	88
3.2.5	Upper Body Radiation Shield, Option	91
3.2.6	Final Work Steps	91
3.3	Part 3: Floor stand	91
3.3.1	Special Requirements	91
3.3.2	Product-specific Safety Information	93
3.3.3	Mechanical Safety	93
3.3.4	Electrical Safety, Safety Devices	98
3.3.5	Maintenance	100
3.3.6	Operating Values, Inspection	103
3.3.7	Function Inspection	104
3.3.8	Final Work Steps	104
3.4	Part 4: Patient Table/UPS	104
3.4.1	Special Requirements	104
3.4.2	Patient Table	106
3.4.3	Option, Uninterruptible Power Supply	114
3.4.4	Final Work Steps	116

4	Partial Maintenance	122
4.1	Special Requirements	122
4.1.1	Required Documents	122
4.1.2	Required Tools, Measuring Equipment and Aids	122
4.1.3	Required Materials	123
4.2	General Visual Check	124
4.2.1	Visual Check per IEC 62353	124
4.3	Image System (AXIS)	125
4.3.1	Product-specific Safety Information	125
4.3.2	Mechanical Safety	125
4.3.3	Electrical Safety	125
4.4	Panoramic Display (option)	126
4.4.1	Product-specific Safety Information	126
4.4.2	Mechanical Safety	126
4.4.3	Electrical Safety	126
4.5	IOB/CCOM (RTC), Input/Output Box, Central Computer	127
4.5.1	Product-specific Safety Information	127
4.5.2	Electrical Safety	127
4.6	MCI	127
4.6.1	Product-specific Safety Information	127
4.6.2	Mechanical Safety	127
4.6.3	Preventive Maintenance	128
4.7	TFT Color Monitor (Examination Room)	128
4.7.1	Product-specific Safety Information	128
4.7.2	Mechanical Safety	128
4.7.3	Inspection and Maintenance	128
4.7.4	Upkeep	130
4.8	21" Display (Control Room)	131
4.8.1	Inspection and Maintenance	131
4.8.2	Upkeep	131
4.9	X-ray Generator	132
4.9.1	Product-specific Safety Information	132
4.9.2	Electrical Safety	133
4.9.3	Preventive Maintenance	134
4.10	Tube Cooling Unit	135
4.10.1	Product-specific Safety Information	135
4.10.2	Preventive Maintenance	136
4.11	DCS	140
4.11.1	Product-specific Safety Information	140
4.11.2	Mechanical Safety	140
4.11.3	Preventive Maintenance	142
4.12	Upper Body Radiation Shield, Option	143
4.12.1	Product-specific Safety Information	143
4.12.2	Mechanical Safety	143
4.13	Stand	143
4.13.1	Product-specific Safety Information	143
4.13.2	Maintenance	144
4.13.3	Electrical Safety, Safety Devices	146
4.14	Patient Table	148
4.14.1	Product-specific Safety Information	148
4.14.2	Mechanical Safety	149
4.14.3	Electrical Safety	149
4.14.4	Periodic Preventive Maintenance	151

4.15	Option, Uninterruptible Power Supply	153
4.15.1	Mechanical Safety	154
4.15.2	Electrical Safety	154
4.15.3	Function Inspection	154
4.15.4	Upkeep	154
4.16	Final Work Steps	155
4.16.1	Image Quality	155
4.16.2	Installing the Cover Panels	156
4.16.3	Performing a System Data Backup	157
4.16.4	Electrical Safety	157
5	Procedures	160
5.1	Protective ground wire test meter	160
5.1.1	Order directly	160
5.2	Cleaning the TFT Displays	160
5.2.1	Cleaning Information	160
6	Changes versus the Previous Version	163
	Index	164

3 Full Maintenance

Type of Maintenance:	Maintenance Interval:	Required time:
Full Maintenance (includes partial maintenance)	Every 24 months	Approx. 17.75 hours

Full maintenance is divided into four partial maintenances. These do not have to be performed in sequence.

Partial Maintenance 1:

Component	Full Maintenance (Work Time/Maintenance Interval)
General Visual Check	0.25 hr. / 24 months
Image system	1 hr. / 24 months
User Interfaces	0.75 hr. / 24 months
Wireless footswitch	0.2 hrs. / 24 months
IOB/CCOM, formerly the RTC	0.5 hr. / 24 months
MCI	0.2 hours/24 months
Examination Room Display	0.5 hr. / 24 months
21" Display	0.5 hr. / 24 months
Final Work Steps	0.75 hr. / 24 months
Total Work Time, Partial Maintenance 1:	Approx. 5 hours/24 months

Partial Maintenance 2:

Component	Full Maintenance (Work Time/Maintenance Interval)
X-ray Generator	2 hrs. / 24 months
Upper Body Protection (optional)	0.25 hr. / 24 months
DCS	0.75 hr. / 24 months
Final Work Steps	0.75 hr. / 24 months
Total Work Time, Partial Maintenance 2:	3.75 hrs. / 24 months

Partial Maintenance 3:

Component	Full Maintenance (Work Time/Maintenance Interval)
Floor Stand	3.5 hours / 24 months
Final Work Steps	0.75 hour / 24 months
Total Work Time, Partial Maintenance 3:	4.25 hours / 24 months

Partial Maintenance 4:

Component	Full Maintenance (Work Time/Maintenance Interval)
Patient Table	3.5 hours / 24 months
Uninterruptible Power Supply (optional)	0.5 hour / 24 months
Final Work Steps	0.75 hour / 24 months
Total Work Time, Partial Maintenance 4:	4.75 hours / 24 months

3.1 Part 1 AXIS Smart/User Interfaces/IOB/Displays**3.1.1 Special Requirements****3.1.1.1 Required Documents**

- Angiography CB-DOC DVD, AXAN-000.801.01....

3.1.1.2 Required Tools, Measuring Equipment and Aids

- Protective ground wire test meter. This test meter must meet the specifications of IEC 61557/EN 61557/VDE 0413/ IEC 62353/, for example the SECUTEST® SIII general-purpose test meter or PROFITEST® 0100S-II. Order directly
- ESD equipment 97 02 606 RE999
- Lint-free cleaning cloths n.a.
- Alcohol n.a.
- Cleaning agent for plastic, glass and painted parts such as Hakapur concentrated cleaner (500 g), 96 60 648 RH999
- Brush to remove dust n.a.
- Anti-static vacuum cleaner n.a.
- Luminous density meter, 77 52 848
such as SMfit ACT (luminous density meter)

All tools, test equipment, and aids, with the exception of the "standard installation tool kit" and the general-purpose test meter, are listed and specified in the Service Parts Catalogue (SPC).

3.1.1.3 Required materials

- 1 lithium battery for the imaging system, type CR2032 30 82 146
- NCU fan/battery module (last replacement > 3 years) 10 96 1600

5 Procedures

5.1 Protective ground wire test meter

5.1.1 Order directly

Info:

- The ground wire test meter (test meter) must meet the specifications of IEC 61557 / EN 61557 / VDE 0413 / IEC 62353,

for example the SECUTEST® SIII general-purpose test meter or PROFITEST® 0100S-II; **order directly from:**

GMC-I Messtechnik GmbH

Südwestpark 15

90449 Nuremberg

Germany

Telephone +49 911 8602-0

Fax +49 911 8602-669

e-mail: info@gossenmetrawatt.com

<http://www.gossenmetrawatt.com>

5.2 Cleaning the TFT Displays

5.2.1 Cleaning Information



Information for cleaning agents/disinfectants apply only to the display to which this document refers and not the entire system.



Information about cleaning agents must be observed. If agents other than those listed are used, the display surface can be irreversibly damaged.

Cleaning agents may not remain on the surface for an extended period.

If this is not observed, the warranty is invalidated.

Equipment care, cleaning and disinfection:

- The front of the display is very sensitive to mechanical damage; absolutely avoid scratches, shocks, etc.
- Promptly remove any drops of liquid that fall on it; extended contact with liquids may cause, for example, discoloration or calcium outlines on the surface.

- If there is dirt on the front of the display, clean it with a microfiber cloth and if necessary with a cleaning agent. Clean housing components only with the recommended cleaning agents.
- If a cleaning agent is sprayed directly on the display surface, make sure that any droplets that run down are wiped off with a microfiber cloth before they reach the edge of the panel.

Active Ingredient Class	Tested Cleaning Agents and Disinfectants	Additional examples
Alcohol	Ethyl alcohol, 96%	Hospiset cloth Microzoid liquid
Aldehyde	Melsitt Cidex	Aldasan 2000 Kohsolin Gigasept FF
Chloride derivative	Terralin	Quartamon Med
Disinfectants	Taski DS5001 (Diverseylever Labs) Morning Mist Surfanios Fraicheur Citron (Anios Labs)	
Guanidine derivatives	Lysoformin	
Quaternary compounds	Incidur spray, full strength	
Commercially available dish-washing liquid	Tempo	Fairy Ultra, Pril, Palmolive
Benzene	Petroleum ether	
Pyridine derivatives	Activ Spray, full strength	
Water	Tap water Distilled water	

Tab. 10: General cleaning agents



Effective time:

The maximum time between applying the cleaning / disinfectant liquid and removal (dry wiping, for example)

If the effective time is exceeded, damage to the surfaces can occur; with the panel, this can be permanent screen damage.

Name	Effective time
Terralin protect; 2%	15 min.
Meliseptol rapid, direct	60 min.
Microbac tissues, direct	30 min.
Incidin Plus, 8% volume	15 min.

Name	Effective time
Melsitt, 10% volume	60 min.
Incidur spray, full strength	60 min.
Perform, 3% weight	60 min.
Lysoformin, 2% volume	60 min.
Activ Spray, full strength	60 min.
Terralin, 0.5%	60 min.
Ethyl alcohol, 96%	90 min.
Commercially available dish washing liquid: Tempo	60 min.
Boiling range - benzene, full strength	Must not be used!
Spray disinfection Disinfectant: Nocolyse Quantity: manufacturer's recommendation	Information from manufacturer

Tab. 11: Directly applied cleaning agents

6 Changes versus the Previous Version

Chapter	Section	Changes
Full Maintenance	Panoramic Display	Title changed from Cockpit to Panoramic Display.
Full Maintenance	MCI	A reminder added to keep F2 switched on before replacement of fan/battery module
Full Maintenance	Checking the Operating Values	Step sequences changed.
Full Maintenance	DCS	Added torque required for DCS installation
Full Maintenance	Floor Stand	Replaced Tube Cover Force Sensor Adjustment with Perform Stand Test
Partial Maintenance	Floor Stand	Replaced Tube Cover Force Sensor Adjustment with Perform Stand Test
System Checkup	DCS	Added a reminder to check the DCS tube covers about buckle fastening and its replacement solution
System Checkup	Final Work Steps	Updated the reminder to check the DCS tube covers about buckle fastening and its replacement solution
General Information	Product-specific Safety information	Service key updated
Full Maintenance	Cockpit (option)	Deleted the reference to the Monitor Replacement Document
Full Maintenance	TFT Color Monitor (Examination Room)	Deleted the reference to the Monitor Replacement Document
Partial Maintenance	TFT Color Monitor (Examination Room)	Deleted the reference to the Monitor Replacement Document
Full Maintenance	X-ray generator	UG-related update added for tube adjustment and service switch
Full Maintenance	Tube Cooling Unit	Updated the signal connector for One4all SMC cooling unit
Partial Maintenance	Tube Cooling Unit	Updated the signal connector for One4all SMC cooling unit

Tab. 12: Artis one - Maintenance Instructions

A

Area dose product 79

C

Cooling unit 84, 136

D

Diamentor 79

H

High voltage connections 76, 134

All documents may only be used by authorized personnel for rendering services on Siemens Healthineers Products. Any document in electronic form may be printed once. Copy and distribution of electronic documents and hard copies is prohibited. Offenders will be liable for damages. All other rights are reserved.

siemens-healthineers.com

Siemens Healthineers Headquarters

Siemens Healthineers AG
Siemensstr. 3
91301 Forchheim
Germany
Telephone: +49 9191 18-0
siemens-healthineers.com