## Testosterone Test System Product Code: 3725-300

### 1.0 INTRODUCTION

**Intended Use: The Quantitative Determination of Total Testosterone Concentration in Human Serum or Plasma by a Microplate Enzyme Immunoassay, Colorimetric**

### 2.0 SUMMARY AND EXPLANATION OF THE TEST

Testosterone, (17 $\beta$ -Hydroxy-4-androstene-3-one), a C<sub>19</sub> steroid, is the most potent naturally secreted androgen.<sup>1</sup> In normal post pubertal males, testosterone is secreted primarily by the testes with only a small amount derived from peripheral conversion of 4-Androstene-3, 17-dione (ASD).<sup>2</sup> In adult women, it has been estimated that over 50% of serum testosterone is derived from peripheral conversion of ASD secreted by the adrenal and ovary, with the remainder from direct secretion of testosterone by these glands.

In the male, testosterone is mainly synthesized in the interstitial Leydig cells and the testis, and is regulated by the interstitial cell stimulating hormone (ICSH), or luteinizing hormone (LH) of the anterior pituitary (the female equivalent of ICSH).<sup>3</sup> Testosterone is responsible for the development of secondary sex characteristics, such as the accessory sex organs, the prostate, seminal vesicles and the growth of facial, pubic and auxiliary hair. Testosterone measurements have been very helpful in evaluating hypogonadal states. Increased testosterone levels in males can be found in complete androgen resistance (testicular feminization). Common causes of decreased testosterone levels in males include: hypogonadism, orchidectomy, estrogen therapy, Klinefelter's syndrome, hypopituitarism, and hepatic cirrhosis.<sup>2,4</sup>


In the female, testosterone levels are normally found to be much lower than those encountered in the healthy male. Testosterone in the female comes from three sources. It is secreted in small quantities by both the adrenal glands and the ovaries, and in healthy women 50-60% of the daily testosterone production arises from peripheral metabolism of prohormone, chiefly androstenedione. Common causes of increased serum testosterone levels in females include polycystic ovaries (Stein-Leventhal syndrome), ovarian tumors, adrenal tumors and adrenal hyperplasia. Virilization in women is associated with the administration of androgens and endogenous overproduction of testosterone. There appears to be a correlation between serum testosterone levels and the degree of virilization in women, although approximately 25% of women with varying degrees of virilism have serum testosterone levels that fall within the female reference range.

### 3.0 PRINCIPLE

#### Competitive Enzyme Immunoassay (TYPE 7):


The essential reagents required for an enzyme immunoassay include antibody, enzyme-antigen conjugate and native antigen.

Upon mixing biotinylated antibody, enzyme-antigen conjugate and a serum containing the native antigen, a competition reaction results between the native antigen and the enzyme-antigen conjugate for a limited number of antibody binding sites. The interaction is illustrated by the following equation:

Ab<sub>BtN</sub> = Biotinylated Antibody (Constant Quantity)  
 Ag = Native Antigen (Variable Quantity)  
 EnzAg = Enzyme-antigen Conjugate (Constant Quantity)  
 AgAb<sub>BtN</sub> = Antigen-Antibody Complex  
 EnzAgAb<sub>BtN</sub> = Enzyme-antigen Conjugate -Antibody Complex  
 k<sub>a</sub> = Rate Constant of Association  
 k<sub>-a</sub> = Rate Constant of Disassociation  
 K = k<sub>a</sub> / k<sub>-a</sub> = Equilibrium Constant

A simultaneous reaction between the biotin attached to the antibody and the streptavidin immobilized on the microwell occurs. This effects the separation of the antibody bound fraction after decantation or aspiration.

The enzyme activity in the antibody bound fraction is inversely proportional to the native antigen concentration. By utilizing several different serum references of known antigen concentration, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

### 4.0 REAGENTS

#### Materials Provided:

#### A. Testosterone Calibrators – 1ml/vial – Icons A-G

Seven (7) vials of serum reference for Testosterone at concentrations of 0 (A), 0.1 (B), 0.5 (C), 1.0 (D), 2.5 (E), 5.0 (F) and 12.0 (G) in ng/ml. Store at 2-8°C. A preservative has been added. The calibrators can be expressed in molar concentrations (nM/L) by multiplying by 3.47. For example: 1ng/ml x 3.47 = 3.47 nM/L

#### B. Testosterone Enzyme Reagent – 6.0 ml/vial – Icon E

One (1) ready to use vial of Testosterone (Analog)-horse radish peroxidase (HRP) conjugate in a protein stabilizing matrix with buffer, red dye, preservative, and binding protein inhibitors. Store 2-8°C.

#### C. Testosterone Biotin Reagent – 6.0 ml – Icon V

One (1) vial containing anti-Testosterone biotinylated purified rabbit IgG conjugate in buffer, dye and preservative. Store at 2-8°C.

#### D. Streptavidin Coated Plate – 96 wells – Icon J

One 96-well microplate coated with 1.0 µg/ml streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

#### E. Wash Solution Concentrate – 20ml/vial – Icon L

One (1) vial containing a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.

#### F. Substrate A – 7ml/vial – Icon S<sup>A</sup>

One (1) vial containing tetramethylbenzidine (TMB) in buffer. Store at 2-8°C.

#### G. Substrate B – 7ml/vial – Icon S<sup>B</sup>

One (1) vial containing hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) in buffer. Store at 2-8°C.

#### H. Stop Solution – 8ml/vial – Icon STOP

One (1) vial containing a strong acid (1N HCl). Store at 2-8°C.

#### I. Product Instructions.

**Note 1:** Do not use reagents beyond the kit expiration date.

**Note 2:** Avoid extended exposure to heat and light. **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on the label.**

**Note 3:** Above reagents are for a single 96-well microplate.

#### 4.1 Required But Not Provided:

- Pipette capable of delivering 0.010ml (10µl), 0.050ml(50µl), 0.100ml ( 100µl) volumes with a precision of better than 1.5%.

- Dispenser(s) for repetitive deliveries of 0.50ml (50µl) .0.100ml (100 µl) and 0.350ml (350µl) volumes with a precision of better than 1.5%.
- Microplate washer or a squeeze bottle (optional).
- Microplate Reader with 450nm and 620nm wavelength absorbance capability.
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplate covers for incubation steps.
- Vacuum aspirator (optional) for wash steps.
- Timer.
- Quality control materials.

### 5.0 PRECAUTIONS

#### For In Vitro Diagnostic Use

Not for Internal or External Use in Humans or Animals

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA required tests. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-8395.

**Safe Disposal of kit components must be according to local regulatory and statutory requirement.**

### 6.0 SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood; serum or plasma in type and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain redtop venipuncture tube or (for plasma) in evacuated tube(s) containing heparin. Allow the blood to clot for serum samples. Centrifuge the specimen to separate the serum or plasma from the cells.

**In patients receiving therapy with high biotin doses (i.e. >5mg/day), no sample should be taken until at least 8 hours after the last biotin administration, preferably overnight to ensure fasting sample.**

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.020ml (20µl) of the specimen is required.

### 7.0 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, normal and high range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

### 8.0 REAGENT PREPARATION

- Wash Buffer**  
Dilute contents of wash solution concentrate to 1000ml with distilled or deionized water in a suitable storage container. Diluted buffer can be stored at 2-30°C for up to 60 days.
- Working Substrate Solution - Stable for 1 year.**  
Pour the contents of the amber vial labeled Solution 'A' into the clear vial labeled Solution 'B'. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2 - 8°C.

**Note 1: Do not use the working substrate if it looks blue.**  
**Note 2: Do not use reagents that are contaminated or have bacteria growth.**

### 9.0 TEST PROCEDURE

Before proceeding with the assay, bring all reagents, serum calibrators and controls to room temperature (20 - 27°C).

**\*\*Test Procedure should be performed by a skilled individual or trained professional\*\***

- Format the microplates' wells for each serum reference, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
- Pipette 0.010 ml (10µL) of the appropriate serum reference, control or specimen into the assigned well.
- Add 0.050 ml (50µl) of the ready to use Testosterone Enzyme Reagent to all wells.
- Swirl the microplate gently for 20-30 seconds to mix.
- Add 0.050 ml (50µl) of Testosterone Biotin Reagent to all wells.
- Swirl the microplate gently for 20-30 seconds to mix.
- Cover and incubate for 60 minutes at room temperature.
- Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
- Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100µl) of working substrate solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells.**  
**DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION**
- Incubate at room temperature for fifteen (15) minutes.
- Add 0.050ml (50µl) of stop solution to each well and gently mix for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
- Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. **The results should be read within thirty (30) minutes of adding the stop solution.**

**Note:** Dilute the samples suspected of concentrations higher than 12 ng/ml 1:5 and 1:10 with Testosterone 'O' ng/ml calibrator or female patient sera with a known low value for testosterone.

### 10.0 CALCULATION OF RESULTS

**A dose response curve is used to ascertain the concentration of Testosterone in unknown specimens.**

- Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
- Plot the absorbance for each duplicate serum reference versus the corresponding Testosterone concentration in ng/ml on linear graph paper (do not average the duplicates of the serum references before plotting).
- Connect the points with a best-fit curve.
- To determine the concentration of Testosterone for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (1.764) intersects the dose response curve at (0.57ng/ml) Testosterone concentration (See Figure 1).


**Note:** Computer data reduction software designed for ELISA assays may also be used for the data reduction. **If such software is utilized, the validation of the software should be ascertained.**

**EXAMPLE 1**

Sample I.D.	Well Number	Abs (A)	Mean Abs (B)	Value (ng/ml)
Cal A	A1	2.780	2.787	0
	B1	2.794		
Cal B	C1	2.576	2.611	0.1
	D1	2.646		
Cal C	E1	1.789	1.877	0.5
	F1	1.965		
Cal D	G1	1.391	1.392	1.0
	H1	1.393		
Cal E	A2	0.780	0.788	2.5
	B2	0.796		
Cal F	C2	0.530	0.538	5.0
	D2	0.547		
Cal G	E2	0.301	0.308	12.0
	F2	0.314		
Ctrl 1	G2	1.040	0.760	1.61
	H2	1.045		
Patient	A3	1.751	1.764	0.57
	B3	1.778		

\*The data presented in Example 1 and Figure 1 is for illustration only and **should not** be used in lieu of a standard curve prepared with each assay.

**Figure 1**

**11.0 Q.C. PARAMETERS**

In order for the assay results to be considered valid the following criteria should be met:

- The absorbance (OD) of calibrator 0 ng/ml should be  $\geq 1.3$
- Four out of six quality control pools should be within the established ranges.

**12.0 RISK ANALYSIS**

The MSDS and Risk Analysis Form for this product are available on request from Monobind Inc.

**12.1 Assay Performance**

- It is important that the time of reaction in each well is held constant to achieve reproducible results.
- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential.

Any deviation from Monobind's IFU may yield inaccurate results.

- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.

12. Risk Analysis: as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from [Monobind@monobind.com](mailto:Monobind@monobind.com).

**12.2 Interpretation**

- Measurements and interpretation of results must be performed by a skilled individual or trained professional.** Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- The reagents for the procedure have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays. (*Boscato LM, Stuart MC, 'Heterophilic antibodies: a problem for all immunoassays' Clin. Chem 1988:3427-33*). For diagnostic purposes the results from this assay should be used in combination with clinical examination, patient's history and, all other clinical findings.
- For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability.**
- If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

**13.0 EXPECTED RANGES OF VALUES**

In agreement with established reference intervals<sup>5</sup> for a "normal" adult population, the expected ranges for the Testosterone AccuBind® ELISA Test System are detailed in Table 1.

**TABLE 1**  
Expected Values for Testosterone EIA Test System (ng/ml)

Boys Before Puberty	0.1 – 3.7
Male	2.5 – 10.0
Female	0.2 – 0.95

It is important to keep in mind that establishment of a range of values, which can be expected to be found by a given method for a population of "normal" persons, is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons, each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

**14.0 PERFORMANCE CHARACTERISTICS**

**14.1 Precision**

The within and between assay precision of the Testosterone AccuBind® ELISA Test System were determined by analyses on three different levels of pool control sera. The number, mean values, standard deviation and coefficient of variation for each of these control sera are presented in Table 2 and Table 3.

**TABLE 2**  
Within Assay Precision (Values in ng/ml)

Sample	N	X	$\sigma$	C.V.
Low	22	1.63	0.16	9.8%
Normal	22	9.14	0.44	4.8%
High	22	14.22	0.79	5.6%

**TABLE 3**  
Between Assay Precision (Values in ng/ml)

Sample	N	X	$\sigma$	C.V.
Low	24	1.72	0.16	9.1%
Normal	24	7.06	0.69	9.7%
High	24	13.08	1.03	7.9%

\*As measured in several experiments in duplicate over a ten day period.

**14.2 Sensitivity**

The Testosterone AccuBind® ELISA Test System has a sensitivity of 0.576 pg. This is equivalent to a sample containing a concentration of 0.0576 ng/ml. The sensitivity was ascertained by determining the variability of the 0 ng/ml serum calibrator and using the 2 $\sigma$  (95% certainty) statistic to calculate the minimum dose.

**14.3 Accuracy**

The Testosterone AccuBind® ELISA Test System was compared with a chemiluminescence immunoassay method. Biological specimens from low, normal and high Testosterone level populations were used. The values ranged from 0.29 ng/ml – 21.9ng/ml. The total number of such specimens was 58. The least square regression equation and the correlation coefficient were computed for this Testosterone EIA in comparison with the reference method. The data obtained is displayed in Table 4.

**TABLE 4**  
Least Square Regression Analysis

Method	Mean (x)	Least Square Regression Analysis	Correlation Coefficient
Monobind (y)	3.12	y = -0.265 + 0.944(x)	0.985
Reference (X)	3.02		

Only slight amounts of bias between this method and the reference method are indicated by the closeness of the mean values. The least square regression equation and correlation coefficient indicates excellent method agreement.

**14.4 Specificity**

The % cross reactivity of the testosterone antibody to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of Testosterone needed to displace the same amount of labeled analog.

Substance	Cross Reactivity
Testosterone	1.0000
Androstenedione	0.0009
Dihydrotestosterone	0.0178
Cortisone	<0.0001
Corticosterone	<0.0001
Cortisol	<0.0001
Spirolactone	<0.0001
Progesterone	<0.0001
17 $\alpha$ -OH Progesterone	<0.0001
DHEA sulfate	<0.0001
Estradiol	<0.0001
Estrone	<0.0001
Estril	<0.0001
Hemolysis	<0.0001
Rubella	<0.0001
Lipemia	<0.0001

**15.0 REFERENCES**

- Dorfman, RI and Shipley, RA, ED: Androgens, New York.: John Wiley and Sons, 1956.
- Horton R, Tait JF: Androstenedione production and interconversion rates measured in peripheral blood and studies on the possible site of conversion to testosterone. J.Clin Invest 45: 301-303, 1966.
- Faiman C and Winter, JSD, Reyes, FI, Clin *Obstet Gynaecol*, 3, 467 (1976).
- Sizonenka, PC, *Pediatrician*, 14, 191 (1987).
- Cummings DC, Wall SR: Non sex hormone binding globulin bound testosterone as a marker for hyperandrogenism. J. Clin Endocrinol Metab. 61:873-876, 1985.
- Lashansky, G, et. al., *J Clin Endocrinol Metab*, 58, 674 (1991)
- Tietz, NW, ED: Clinical Guide to Laboratory Tests, 3<sup>rd</sup> ed. Philadelphia, WA Saunders Co, 1995.

Size	96(A)	192(B)
Reagent (fill)	A) 1ml set	1ml set
	B) 1 (6ml)	2 (6ml)
	C) 1 (6ml)	2 (6ml)
	D) 1 plate	2 plates
	E) 1 (20ml)	1 (20ml)
	F) 1 (7ml)	2 (7ml)
	G) 1 (7ml)	2 (7ml)
	H) 1 (8ml)	2 (8ml)

For Orders and Inquires, please contact


Tel: +1 949.951.2665 Mail: info@monobind.com  
Fax: +1 949.951.3539 Fax: www.monobind.com

Please visit our website to learn more about our products and services.

**Glossary of Symbols**  
(EN 980/ISO 15223)

## Thyroglobulin (Tg) Test System Product Code: 2225-300

### 1.0 INTRODUCTION

**Intended Use: The Quantitative Determination of Thyroglobulin (Tg) Concentration in Human Serum by a Microplate Enzyme immunoassay, Colorimetric**

### 2.0 SUMMARY AND EXPLANATION OF THE TEST

Human thyroglobulin (Tg) is a large glycoprotein (660 kD) that is stored in the follicular colloid of the thyroid gland. It functions as a prohormone in the intrathyroid synthesis of primary thyroid hormones like *Triiodothyronine (T3)* and *Thyroxine (T4)*.

Tg is elevated in thyroid follicular and papillary carcinoma, thyroid adenoma, subacute thyroiditis, Hashimoto's thyroiditis and Graves Disease. Tg levels are found to be normal in patients with medullary thyroid carcinoma. Serial measurements of Tg is most useful in detecting recurrence of differentiated thyroid carcinoma following surgical resection or radioactive iodine ablation. Tg determination is used as an adjunct to iodine scanning but not as a replacement for it. Assessment of Tg levels aids in management of infants with congenital hypothyroidism.


Tg determination has been done with various methods using direct competitive binding RIA and double antibody sandwich IRMA or ELISA, of which latter is more useful. All these methods suffer from interference by endogenous autoantibodies to Tg. It is useful to determine the effect of autoantibodies before screening such patients for levels of Tg. *Monobind* provides Tg autoantibody ELISA to rule out such interference. (**Please see *Monobind Anti-Tg AccuBind® ELISA Test System, Product Code: 1025-300***).

### 3.0 PRINCIPLE

#### Immunoenzymometric sequential assay (TYPE 4):

The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitope recognition, in excess, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on the well and exogenously added biotinylated monoclonal Thyroglobulin antibody.

When monoclonal biotinylated antibody is mixed with a serum containing the Tg antigen, a reaction results between the Tg antigen and the antibody, to form an antibody-antigen complex. Simultaneously the biotin attached to the antibody binds to the streptavidin coated on the microwells resulting in immobilization of the complex. The interaction is illustrated by the following equation:


$B^{in}Ab_{(m)}$  = Biotinylated Monoclonal Antibody (Excess Quantity)  
 $Ag_{(Tg)}$  = Native Antigen (Variable Quantity)  
 $Ag_{(Tg)} - B^{in}Ab_{(m)}$  = Antigen-Antibody complex (Variable Quan.)

$k_a$  = Rate Constant of Association  
 $k_{-a}$  = Rate Constant of Disassociation

$Ag_{(Tg)} - B^{in}Ab_{(m)} + Streptavidin_{C.W.} \Rightarrow Immobilized\ complex\ (IC)$   
 $Streptavidin_{C.W.}$  = Streptavidin immobilized on well  
 $Immobilized\ complex\ (IC)$  = Ag-Ab bound to the well

After a suitable incubation period, the antibody-antigen bound fraction is separated from unbound antigen by decantation or aspiration. Another antibody (directed at a different epitope) labeled with an enzyme is added. Another interaction occurs to form an enzyme labeled antibody-antigen-biotinylated-antibody complex on the surface of the wells. Excess enzyme is washed off via a wash step. A suitable substrate is added to produce color measurable with the use of a microplate spectrophotometer. The enzyme activity on the well is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen concentration, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

$EnzAb_{(Tg)}$  = Enzyme labeled Antibody (Excess Quantity)  
 $EnzAb_{(Tg)} - IC$  = Antigen-Antibodies Complex  
 $k_b$  = Rate Constant of Association  
 $k_{-b}$  = Rate Constant of Disassociation

### 4.0 REAGENTS

#### Materials Provided:

#### A. Thyroglobulin Calibrators – 1.0 ml/vial – Icons A - F

Six (6) vials of references for Thyroglobulin antigen (Tg) at levels of **0(A), 2.0 (B), 10.0(C), 40(D), 100(E), and 250(F) ng/ml**. A preservative has been added.

**Note:** There is no known, internationally accepted thyroglobulin standard available. The Tg used in the serum based calibrators is a highly purified (98+% pure) human Tg preparation that is calibrated gravimetrically against the reference material obtained from Community Bureau of Reference # CRM 457.

#### B. x-Tg Biotin Reagent – 13ml/vial – Icon ▽

One (1) vial contains biotinylated anti-Tg monoclonal mouse IgG in buffer, dye, and preservative. Store at 2-8°C.

#### C. Tg Enzyme Reagent – 13 ml/vial - Icon ⊕

One (1) vial contains anti-thyroglobulin IgG labeled with horseradish peroxidase (HRP) in buffer, dye, and preservative. Store at 2-8°C.

#### D. Streptavidin Coated Plate – 96 wells – Icon ↓

One 96-well (creat well modules) microplate coated with streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

#### E. Wash Solution Concentrate – 20 ml/vial – Icon ♣

One (1) vial contains surfactant in buffered saline. A preservative has been added. Store at 2-8°C.

#### F. Substrate Reagent – 12ml/vial – Icon S<sup>N</sup>

One (1) vial contains tetramethylbenzidine (TMB) and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) in buffer. Store at 2-8°C.

#### G. Stop Solution – 8ml/vial – Icon ⊖

One (1) vial contains a strong acid (0.5M H<sub>2</sub>SO<sub>4</sub>). Store at 2-8°C.

#### H. Product Instructions.

**Note 1:** Do not use reagents beyond the kit expiration date.

**Note 2:** Avoid extended exposure to heat and light. **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on the label.**

**Note 3:** Above reagents are for a single 96-well microplate.

#### 4.1 Required But Not Provided:

- Pipette(s) capable of delivering 0.050ml (50µl) and 0.100ml (100µl) volumes with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.100ml (100µl) and 0.350ml (350µl) volumes with a precision of better than 1.5%
- Microplate washer or a squeeze bottle (optional).
- Microplate Reader with 450nm and 620nm wavelength absorbance capability
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplate cover for incubation steps.
- Vacuum aspirator (optional) for wash steps.

- Timer.
- Storage container for storage of wash buffer.
- Distilled or deionized water.
- Quality Control Materials.

### 5.0 PRECAUTIONS

**For In Vitro Diagnostic Use  
Not for Internal or External Use in Humans or Animals**

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface antigen, HIV 1&2 and HCV antibodies by FDA required tests. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS.

**Safe Disposal of kit components must be according to local regulatory and statutory requirement.**

### 6.0 SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood serum in type, and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain red-top venipuncture tube without additives or gel barrier. Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells.

**In patients receiving therapy with high biotin doses (i.e. >5mg/day), no sample should be taken until at least 8 hours after the last biotin administration, preferably overnight to ensure fasting sample.**

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.100ml (100µl) of the specimen is required.

### 7.0 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, normal and elevated range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

### 8.0 REAGENT PREPARATION

#### 1. Wash Buffer

Dilute contents of wash concentrate to 1000ml with distilled or deionized water in a suitable storage container Diluted buffer can be stored at room temperature (2-30°C) for up to 60 days.

**Note: Do not use reagents that are contaminated or have bacteria growth.**

#### 9.0 TEST PROCEDURE (Time 4hr 15min)

*Before proceeding with the assay, bring all reagents, serum reference calibrator and controls to room temperature (20-27°C).*

**\*\*Test Procedure should be performed by a skilled individual or trained professional\*\***

- Format the microplates' wells for each calibrator, control and patient sample to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
- Pipette 0.050 ml (50µl) of the appropriate calibrators, controls and samples into the assigned wells.
- Add 0.100 ml (100µl) of the x-Tg Biotin Reagent to each well. **It is very important to dispense all reagents close to the bottom of the microwell.**

- Swirl the microplate gently for 20-30 seconds to mix. Cover with a plastic wrap or microplate cover.
- Incubate for 2 hours at room temperature.
- Discard the contents of the microplate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
- Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is used, fill each well to the top by squeezing the container (Avoiding air bubbles). Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100µl) of Tg Enzyme Reagent to all wells **DO NOT SHAKE THE PLATE AFTER ENZYME ADDITION**
- Cover with a plastic wrap. Incubate at room temperature for 120 minutes.
- Discard the contents of the microplate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
- Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is used, fill each well to the top by squeezing the container (Avoiding air bubbles). Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100 µl) of substrate to all wells. **Always add reagents in the same order to minimize reaction time differences between wells.**  
**DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION**
- Cover with a plastic wrap or microplate cover. Incubate at room temperature for 15 minutes.
- Add 0.050ml (50µl) of stop solution to each well and mix gently for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
- Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. **The results should be read within fifteen (15) minutes of adding the stop solution.**

#### 9.1 ALTERNATE PROCEDURE (Time 2hr 15min)

**This procedure can be used with the help of a laboratory hematology shaker.**

- Format the microplates' wells for each calibrator, control and patient sample to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
- Pipette 0.050 ml (50µl) of the appropriate calibrators, controls and samples into the assigned wells.
- Add 0.100 ml (100µl) of the biotin labeled monoclonal antibody to each well. **It is very important to dispense all reagents close to the bottom of the microwell and swirl to mix.**
- Incubate at room temperature for 1 hour while shaking constantly on a hematology shaker at 150 RPM.
- Discard the contents of the microplate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
- Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is used, fill each well to the top by squeezing the container (Avoiding air bubbles). Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100µl) of Tg Enzyme Reagent to all wells
- Incubate at room temperature for 1 hour while shaking constantly on a hematology shaker at 150 RPM.
- Repeat steps 5-6 as described in the "Test Procedure" above.
- Follow steps 11-14 to develop color and measure.

### 10.0 CALCULATION OF RESULTS

**A dose response curve is used to ascertain the concentration of human thyroglobulin (Tg) in unknown specimens.**

- Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
- Plot the absorbance for each duplicate serum reference versus the corresponding Tg concentration in ng/ml on linear graph

paper (do not average the duplicates of the serum references before plotting).

- Draw the best-fit curve through the plotted points.
- To determine the concentration of Tg for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (0.424) intersects the dose response curve at 25.2 ng/ml Tg concentration (See Figure 1).


**Note:** Computer data reduction software designed for ELISA assays may also be used for the data reduction. **If such software is utilized, the validation of the software should be ascertained.**

#### EXAMPLE 1

(Data is based on 9.0 Test Procedure - time 4 hr. 15min)

Sample I.D.	Well Number	Abs (A)	Mean Abs (B)	Value (ng/ml)
Cal A	A1	0.047	0.047	0
	B1	0.047		
Cal B	C1	0.093	0.091	2
	D1	0.090		
Cal C	E1	0.221	0.217	10
	F1	0.214		
Cal D	G1	0.612	0.625	40
	H1	0.634		
Cal E	A2	1.343	1.339	100
	B2	1.335		
Cal F	C2	2.596	2.577	250
	D2	2.557		
Cont 1	E2	0.142	0.146	4.99
	F2	0.150		
Cont 2	G2	1.622	0.876	125.0
	H2	1.566		
Patient 1	A3	0.426	0.424	25.2
	B3	0.422		

Figure 1

\*The data presented in Example 1 and Figure 1 are for illustration only and **should not** be used in lieu of a dose response curve prepared with each assay.

#### 11.0 Q. C. PARAMETERS

In order for the assay results to be considered valid the following criteria should be met.

- The absorbance (OD) of calibrator A should be  $\leq 0.06$ .
- The absorbance (OD) of calibrator F should be  $\geq 1.3$ .
- Four out of six quality control pools should be within the established ranges.

#### 12.0 RISK ANALYSIS

The MSDS and Risk Analysis Form for this product are available on request from Monobind Inc.

#### 12.1 Assay Performance

- It is important that the time of reaction in each well is held constant to achieve reproducible results.

- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. The substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Patient samples with thyroglobulin concentrations above 250 ng/ml may be diluted with the zero calibrator and re-assayed. Multiply the value obtained by the dilution factor to obtain the corrected value.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind IFU may yield inaccurate results.
- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
- Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from [Monobind@monobind.com](mailto:Monobind@monobind.com).

#### 12.2 Interpretation

- Measurements and interpretation of results must be performed by a skilled individual or trained professional.
- Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- The reagents for the test system have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of Immunoassays (Boscato LM, Stuart MC "Heterophilic antibodies: a problem for all immunoassays" *Clin. Chem* (1988:3427-33). For diagnostic purposes, the results from this assay should be in combination with clinical examination, patient history, and all other clinical findings.
- For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability.**
- If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

#### 13.0 EXPECTED RANGE OF VALUES

Based on the clinical data gathered by Monobind in concordance with the published literature a normal range was established.

Table 1: Expected Values for Tg	
POPULATION	RANGE
Adult	3.5 – 56 ng/ml

Tg is found to be elevated in patients with thyroid follicular and papillary carcinoma, thyroid adenoma, subacute thyroiditis, Hashimoto's thyroiditis and Graves' disease. Low levels of Tg are an indication of thyrotoxicosis factitia.

It is important to keep in mind that any normal range establishment is dependent upon a multiplicity of factors like the specificity of the method, the locale, the population tested and the precision of the method in the hands of technicians. For these reasons each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the technicians using the

method with a population indigenous to the area in which the laboratory is located.

#### 14.0 PERFORMANCE CHARACTERISTICS

##### 14.1 Precision

The within and between assay precisions of the Thyroglobulin AccuBind® ELISA test system were determined by analyses on three different levels of pool control sera. The number (N), mean value (X), standard deviation ( $\sigma$ ) and coefficient of variation (C.V.) for each of these control sera are presented in Tables 2 and 3.

TABLE 2				
Within Assay Precision (Values in ng/ml)				
Sample	N	X	$\sigma$	C.V.
Pool 1	22	6.2	0.41	6.6%
Pool 2	22	64.4	2.23	3.6%
Pool 3	22	194.1	8.17	4.2%

TABLE 3				
Between Assay Precision* (Values in ng/ml)				
Sample	N	X	$\sigma$	C.V.
Pool 1	10	5.8	0.52	9.0%
Pool 2	10	62.2	3.82	6.1%
Pool 3	10	192.3	10.90	5.7%

\*As measured in ten experiments in duplicate over seven days.

##### 14.2 Sensitivity

The analytical sensitivity (detection limit) was ascertained by determining the variability of the 0 ng/ml serum calibrator and using the  $2\sigma$  (95% certainty) statistic to calculate the minimum dose. The assay sensitivity was found to be 0.44 ng/ml.

##### 14.3 Accuracy

The Tg AccuBind® ELISA test system was compared with a reference coated tube radioimmunoassay (IRMA) assay. Biological specimens from symptomatic and asymptomatic populations were used. The data obtained is displayed in Table 4.

TABLE 4			
Method	Mean (x)	Least Square Regression Analysis	Correlation Coefficient
Monobind	13.6	$y = 2.55 + 0.908(x)$	0.975
Reference	11.4		

The least square regression equation and correlation coefficient indicates excellent method agreement.

##### 14.4 Specificity

The cross-reactivity of the Tg AccuBind® ELISA method to selected substances was evaluated by adding the interfering substance(s) to a serum matrix at the following concentration(s). The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of Thyroglobulin needed to produce the same absorbance.

Substance	Concentration	Cross Reactivity
Thyroglobulin	100 ng/ml	100.0%
Triiodothyronine	1000 ng/dl	N/D
Thyroxine	1000 ng/ml	N/D
TBG	100 ng/ml	N/D

##### 14.5 High Dose Effect

Since the assay is sequential in design, high concentrations of Tg do not show the hook effect. Samples with concentrations over 50,000 ng/ml demonstrated extremely high levels of absorbance.

#### 15.0 REFERENCES

- Beever K, Bradbury J, Phillips D, et al, "Highly sensitive assays of autoantibodies to Thyroglobulin and Thyroid Peroxidase", *Clin Chem*, 35, 1949-1954 (1989).
- Ladenson PW, "Optimal laboratory testing for diagnosis and monitoring of thyroid nodules, goiter, and thyroid cancer", *Clin Chem*, 42, 183-187 (1996).
- Mayo Medical Laboratories: test Catalog, Rochester, MN (1997).
- Spencer CA, TakeUCHO M, Kazarosyn M, "Current status and performance goals for serum thyroglobulin assays", *Clin Chem*, 42, 164-173 (1996).
- Tietz N. Ed: *Clinical Guide to Laboratory Tests*. 3<sup>rd</sup> Ed. Philadelphia. Saunders (1995).

- Surks, MI, Chopra, IJ, Mariash, CN, "American Thyroid Association guidelines for use of laboratory tests in thyroid disorders", *JAMA*, 263, 1529-1532 (1990).
- Ng, M., Rajna, A, Khalid, B, "Enzyme immunoassay for simultaneous measurements of autoantibodies against thyroglobulin and thyroid microsomes in serum", *Clin Chem*, 33, 2286-2288 (1987).
- Spencer CA, TakeUCHO M, Kazarosyn M, Wang CC, Guttler RB, Singer PA, et al, "Serum thyroglobulin autoantibodies; prevalence, influence on serum thyroglobulin measurements, and prognostic significance in patients with differentiated thyroid carcinoma", *J Clin Endocrinol Metab*, 83, 1121-127 (1998).
- Spencer CA, LoPresti JS, Fatemi S, Nicoloff JT, "Detection of residual and recurrent differentiated thyroid carcinoma by serum thyroglobulin measurements", *Thyroid*, 9, 435-41 (1999).
- Schlumberger M, Baudin E, "Serum thyroglobulin determinations in the follow up of patients with differentiated thyroid carcinoma", *Eur J. Endocrinol*, 138, 249-252 (1998).

Effective Date: 2019-Jul-16 Rev. 4  
MP2225


DCO: 1353  
Product Code: 2225-350

Size	96(A)	192(B)	
Reagent (fill)	A)	1ml set	1ml set
	B)	1 (13ml)	2 (13ml)
	C)	1 (13ml)	2 (13ml)
	D)	1 plate	2 plates
	E)	1 (20ml)	1 (20ml)
	F)	1 (12ml)	2 (12ml)
	G)	1 (8ml)	2 (8ml)

For Orders and Inquires, please contact

**Monobind Inc.**  
100 North Pointe Drive  
Lake Forest, CA 92630 USA


Tel: +1 949.951.2655 Mail: [info@monobind.com](mailto:info@monobind.com)  
Fax: +1 949.951.3539 Fax: [www.monobind.com](http://www.monobind.com)

CEpartner4U, Esdoornlaan 13  
3951 DBMaarn, The Netherlands  
[www.cepartner4u.eu](http://www.cepartner4u.eu)

Please visit our website to learn more about our products and services.

#### Glossary of Symbols (EN 980/ISO 15223)

## Thyrotropin (TSH) Test System Product Code: 325-300

### 1.0 INTRODUCTION

**Intended Use: The Quantitative Determination of Thyrotropin Concentration in Human Serum by a Microplate Enzyme Immunoassay, Colorimetric**

### 2.0 SUMMARY AND EXPLANATION OF THE TEST

Measurement of the serum concentration of thyrotropin (TSH), a glycoprotein with a molecular weight of 28,000 Daltons and secreted from the anterior pituitary, is generally regarded as the most sensitive indicator available for the diagnosis of primary and secondary (pituitary) hypothyroidism.<sup>1,2</sup> The structure of human TSH is similar to that of the pituitary and placental gonadotropins, consisting of an 89-amino acid  $\alpha$ -subunit which is similar or identical between these hormones and a 115-amino acid  $\beta$ -subunit, which apparently confers hormonal specificity. The production of the 2 subunits is separately regulated with apparent excess production of the  $\alpha$ -subunit. The TSH molecule has a linear structure consisting of the protein core with carbohydrate side chains; the latter accounts for 16% of the molecular weight.

TSH measurements are equally useful in differentiating secondary and tertiary (hypothalamic) hypothyroidism from the primary thyroid disease. TSH release from the pituitary is regulated by thyrotropin releasing factor (TRH), which is secreted by the hypothalamus, and by direct action of T4 and triiodothyronine (T3), the thyroid hormones, at the pituitary. Increase levels of T3 and T4 reduces the response of the pituitary to the stimulatory effects of TRH. In secondary and tertiary hypothyroidism, concentrations of T4 are usually low and TSH levels are generally low or normal. Either pituitary TSH deficiency (secondary hypothyroidism) or insufficiency of stimulation of the pituitary by TRH (tertiary hypothyroidism) causes this. The TRH stimulation test differentiates these conditions. In secondary hypothyroidism, TSH response to TRH is blunted while a normal or delayed response is obtained in tertiary hypothyroidism.

Further, the advent of immunoassay methods has provided the laboratory with sufficient sensitivity to enable the differentiating of hyperthyroidism from euthyroid population and extending the usefulness of TSH measurements. This method is a second-generation assay, which provides the means for discrimination in the hyperthyroid-euthyroid range. The functional sensitivity (<20% between assay CV) of the one-hour procedure is 0.195  $\mu$ U/ml while the two-hour procedure has a functional sensitivity of 0.095  $\mu$ U/ml.<sup>3</sup>

In this method, TSH calibrator, patient specimen or control is first added to a streptavidin coated well. Biotinylated monoclonal and enzyme labeled antibodies are added and the reactants mixed. Reaction between the various TSH antibodies and native TSH forms a sandwich complex that binds with the streptavidin coated to the well.

After the completion of the required incubation period, the antibody bound enzyme-thyrotropin conjugate is separated from

the unbound enzyme-thyrotropin conjugate by aspiration or decantation. The activity of the enzyme present on the surface of the well is quantitated by reaction with a suitable substrate to produce color.


The employment of several serum references of known thyrotropin levels permits construction of a dose response curve of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with thyrotropin concentration.

### 3.0 PRINCIPLE

#### Immunoassay (TYPE 3):

The essential reagents required for an immunoassay include high affinity and specificity antibodies (enzyme conjugated and immobilized), with different and distinct epitope recognition, in excess, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on the well and exogenously added biotinylated monoclonal anti-TSH antibody.

Upon mixing monoclonal biotinylated antibody, the enzyme-labeled antibody and a serum containing the native antigen, reaction results between the native antigen and the antibodies, without competition or steric hindrance, to form a soluble sandwich complex. The interaction is illustrated by the following equation:

$\text{BiotAb}_{(m)}$  = Biotinylated Monoclonal Antibody (Excess Quantity)  
 $\text{Ag}_{\text{TSH}}$  = Native Antigen (Variable Quantity)  
 $\text{EnzAb}_{(p)}$  = Enzyme-Polyclonal Antibody (Excess Quantity)  
 $\text{EnzAb}_{(p)} - \text{Ag}_{\text{TSH}} - \text{BiotAb}_{(m)}$  = Antigen-Antibodies Sandwich Complex  
 $k_a$  = Rate Constant of Association  
 $k_a$  = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:  
 $\text{EnzAb}_{(p)} - \text{Ag}_{\text{TSH}} - \text{BiotAb}_{(m)} + \text{Streptavidin}_{\text{CW}} \rightarrow \text{immobilized complex}$   
 $\text{Streptavidin}_{\text{CW}}$  = Streptavidin immobilized on well  
 Immobilized complex = sandwich complex bound to the well surface

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

### 4.0 REAGENTS

#### Materials Provided:

##### A. TSH Calibrators – 1ml/vial - Icons A-G

Seven (7) vials of references for TSH Antigen at levels of 0(A), 0.5(B), 2.5(C), 5.0(D), 10(E), 20(F) and 40(G)  $\mu$ U/ml. Store at 2-8°C. A preservative has been added.

**Note:** The calibrators, human serum based, were calibrated using a reference preparation, which was assayed against the WHO 2<sup>nd</sup> IRP 80/558.

##### B. TSH Enzyme Reagent – 13ml/vial - Icon $\text{\textcircled{E}}$

One (1) vial containing enzyme labeled affinity purified polyclonal goat antibody, biotinylated monoclonal mouse IgG in buffer, dye, and preservative. Store at 2-8°C.

##### C. Streptavidin Coated Plate – 96 wells - Icon $\text{\textcircled{J}}$

One 96-well microplate coated with streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

##### D. Wash Solution Concentrate – 20 ml/ml - Icon $\text{\textcircled{L}}$

One (1) vial containing a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.

##### E. Substrate A – 7ml/vial - Icon $\text{\textcircled{S}}$

One (1) vial containing tetramethylbenzidine (TMB) in buffer. Store at 2-8°C.

##### F. Substrate B – 7ml/vial - Icon $\text{\textcircled{B}}$

One (1) vial containing hydrogen peroxide ( $\text{H}_2\text{O}_2$ ) in buffer. Store at 2-8°C.

### G. Stop Solution – 8ml/vial - Icon $\text{\textcircled{STOP}}$

One (1) vial containing a strong acid (1N HCl). Store at 2-8°C.

### H. Product Instructions.

**Note 1:** Do not use reagents beyond the kit expiration date.

**Note 2:** Avoid extended exposure to heat and light. **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on the label.**

**Note 3:** Above reagents are for a single 96-well microplate.

#### 4.1 Required But Not Provided:

- Pipette(s) capable of delivering 0.050ml (50 $\mu$ l) and 0.100ml (100 $\mu$ l) volumes with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.100ml (100 $\mu$ l) and 0.350ml (350 $\mu$ l) volumes with a precision of better than 1.5% (optional).
- Microplate washer or a squeeze bottle (optional).
- Microplate Reader with 450nm and 620nm wavelength absorbance capability.
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplate cover for incubation steps.
- Vacuum aspirator (optional) for wash steps.
- Timer.
- Storage container for storage of wash buffer.
- Distilled or deionized water.
- Quality Control Materials.

### 5.0 PRECAUTIONS

**For In Vitro Diagnostic Use  
Not for Internal or External Use in Humans or Animals**

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface antigen, HIV 1&2 and HCV antibodies by FDA required tests. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS.

**Safe disposal of kit components must be according to local regulatory and statutory requirement.**

### 6.0 SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood, serum in type, and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain redtop venipuncture tube without additives or gel barrier. Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells.

**In patients receiving therapy with high biotin doses (i.e. >5mg/day), no sample should be taken until at least 8 hours after the last biotin administration, preferably overnight to ensure fasting sample.**

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, (100 $\mu$ l) 0.100 ml of the specimen is required.

### 7.0 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, normal, and elevated range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the dose response curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in

experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

### 8.0 REAGENT PREPARATION

#### 1. Wash Buffer

Dilute contents of wash concentrate to 1000ml with distilled or de-ionized water in a suitable storage container. Store at 2-30°C for up to 60 days.

#### 2. Working Substrate Solution – Stable for one year

Pour the contents of the amber vial labeled Solution 'A' into the clear vial labeled Solution 'B'. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2 - 8°C.

**Note1: Do not use the working substrate if it looks blue.**

**Note 2: Do not use reagents that are contaminated or have bacteria growth.**

### 9.0 TEST PROCEDURE

*Before proceeding with the assay, bring all reagents, serum reference calibrators and controls to room temperature (20-27°C).*

**\*\*Test Procedure should be performed by a skilled individual or trained professional\*\***

- Format the microplates' wells for each serum reference calibrator, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
- Pipette 0.050 ml (50 $\mu$ l) of the appropriate serum reference, control or specimen into the assigned well.
- Add 0.100 ml (100 $\mu$ l) of the TSH Enzyme Reagent to each well. **It is very important to dispense all reagents close to the bottom of the coated well.**
- Swirl the microplate gently for 20-30 seconds to mix and cover.
- Incubate 60 minutes at room temperature. \*\*
- Discard the contents of the microplate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
- Add 0.350ml (350 $\mu$ l) of wash buffer (see Reagent Preparation Section) decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100 $\mu$ l) of working substrate solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells.**
- DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION**
- Incubate at room temperature for fifteen (15) minutes.
- Add 0.050ml (50 $\mu$ l) of stop solution to each well and mix gently for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
- Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. **The results should be read within thirty (30) minutes of adding the stop solution.**

\*\* For better low-end sensitivity (< 0.5 $\mu$ U/ml), incubate 120 minutes at room temperature. The 40 $\mu$ U/ml calibrator should be excluded since absorbance over 3.0 units will be experienced. Follow the remaining steps.

**Note:** Dilute samples reading over 40  $\mu$ U/ml by 1:5 and 1:10 with TSH '0' Calibrator. Multiply the results by the dilution factor to obtain accurate results.

### 10.0 CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of thyrotropin in unknown specimens.

- Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1
- Plot the absorbance for each duplicate serum reference versus the corresponding TSH concentration in  $\mu$ U/ml on linear graph paper (do not average the duplicates of the serum references before plotting).


- Draw the best-fit curve through the plotted points.
- To determine the concentration of TSH for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in  $\mu\text{IU/ml}$ ) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (0.775) intersects the dose response curve at (7.66  $\mu\text{IU/ml}$ ) TSH concentration (See Figure 1).

**Note:** Computer data reduction software designed for ELISA assay may also be used for the data reduction. If such software is utilized, the validation of the software should be ascertained.

**EXAMPLE 1**

Sample I.D.	Well Number	Abs	Mean Abs	Value ( $\mu\text{IU/ml}$ )
Cal A	A1	0.018	0.019	0
	B1	0.021		
Cal B	C1	0.076	0.079	0.5
	D1	0.082		
Cal C	E1	0.302	0.298	2.5
	F1	0.293		
Cal D	G1	0.556	0.567	5.0
	H1	0.577		
Cal E	A2	0.926	0.921	10
	B2	0.916		
Cal F	C2	1.610	1.619	20
	D2	1.629		
Cal G	E2	2.694	2.671	40
	F2	2.647		
Control	G2	0.800	0.775	7.66
	H2	0.751		
Patient	A3	1.391	1.383	16.65
	B3	1.375		

**Figure 1**

\*The data presented in Example 1 and Figure 1 are for illustration only and **should not** be used in lieu of a dose response curve prepared with each assay.

### 11.0 Q.C. PARAMETERS

In order for the assay results to be considered valid the following criteria should be met:

- The absorbance of calibrator 'G' (40  $\mu\text{IU/ml}$ ) should be  $\geq 1.3$ .
- Four out of six quality control pools should be within the established ranges.

### 12.0 RISK ANALYSIS

The MSDS and Risk Analysis Form for this product is available on request from Monobind Inc.

#### 12.1 Assay Performance

- It is important that the time of reaction in each well is held constant to achieve reproducible results.
- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in

the same sequence to eliminate any time-deviation during reaction.

- Plate readers measure vertically. Do not touch the bottom of the wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind IFU may yield inaccurate results.
- Patient specimens with TSH concentrations over 40  $\mu\text{IU/ml}$  may be diluted (1:5 or 1:10) with the '0' calibrator and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor.
- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
- Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from [Monobind@monobind.com](mailto:Monobind@monobind.com).

#### 12.2 Interpretation

**1. Measurement and interpretation of results must be performed by a skilled individual or trained professional.**

- Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- The reagents for the test system have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays (Boscato LM, Stuart MC. 'Heterophilic antibodies: a problem for all immunoassays' Clin. Chem. 1988:3427-33). For diagnostic purposes, the results from this assay should be in combination with clinical examination, patient history and all other clinical findings. For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability.**
- If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.
- Serum TSH concentration is dependent upon a multiplicity of factors: hypothalamus gland function, thyroid gland function, and the responsiveness of pituitary to TRH. **Thus, thyrotropin concentration alone is not sufficient to assess clinical status.**

- Serum TSH values may be elevated by pharmacological intervention. Domperidone, amiodazon, iodide, phenobarbital, and phenytoin have been reported to increase TSH levels.
- A decrease in thyrotropin values has been reported with the administration of propranolol, methimazol, dopamine and d-thyroxine.<sup>4</sup>
- Genetic variations or degradation of intact TSH into subunits may affect the binding characteristics of the antibodies and influence the final result. Such samples normally exhibit different results among various assay systems due to the reactivity of the antibodies involved.

**"NOT INTENDED FOR NEWBORN SCREENING"**

#### 13.0 EXPECTED RANGES OF VALUES

A study of euthyroid adult population was undertaken to determine expected values for the TSH AccuBind® ELISA Test System. The number and determined range are given in Table 1. A nonparametric method (95% Percentile Estimate) was used.

TABLE 1 Expected Values for the TSH ELISA Test System (in $\mu\text{IU/ml}$ )		
Number	139	2.5 Percentile-70% Conf Int
Low Normal	0.39	Low Range 0.28 - 0.53
High Normal	6.16	High Range 5.60 - 6.82

It is important to keep in mind that establishment of a range of values which can be expected to be found by a given method for a population of "normal"-persons is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

#### 14.0 PERFORMANCE CHARACTERISTICS

##### 14.1 Precision

The within and between assay precisions of the TSH AccuBind® test system were determined by analyses on three different levels of pool control sera. The number (N), mean (X) value, standard deviation ( $\sigma$ ) and coefficient of variation (C.V.) for each of these control sera are presented in Table 2 and Table 3.

**TABLE 2**  
Within Assay Precision (Values in  $\mu\text{IU/ml}$ )

Sample	N	X	$\sigma$	C.V.
Pool 1	24	0.37	0.03	8.1%
Pool 2	24	6.75	0.43	6.4%
Pool 3	24	29.30	1.94	6.6%

**TABLE 3**  
Between Assay Precision\* (Values in  $\mu\text{IU/ml}$ )

Sample	N	X	$\sigma$	C.V.
Pool 1	10	0.43	0.04	9.3%
Pool 2	10	6.80	0.54	7.9%
Pool 3	10	28.40	1.67	5.9%

\*As measured in ten experiments in duplicate over seven days.

##### 14.2 Sensitivity

The sensitivity (detection limit) was ascertained by determining the variability of the 0  $\mu\text{IU/ml}$  serum calibrator and using the 2 $\sigma$  (95% certainty) statistic to calculate the minimum dose:

For 1 hr incubation = 0.078  $\mu\text{IU/ml}$   
For 2 hr incubation = 0.027  $\mu\text{IU/ml}$

##### 14.3 Accuracy

The TSH AccuBind® ELISA test system was compared with a reference immunochemiluminescence assay. Biological specimens from hypothyroid, euthyroid and hyperthyroid populations were used (The values ranged from 0.01  $\mu\text{IU/ml}$  - 61  $\mu\text{IU/ml}$ ). The total number of such specimens was 241. The least square regression equation and the correlation coefficient were computed for the TSH AccuBind® ELISA method in comparison with the reference method. The data obtained is displayed in Table 4.

**TABLE 4**

Method	Mean (x)	Least Square Regression Analysis	Correlation Coefficient
Monobind	4.54	$y = 0.47 + 0.968 (x)$	0.995
Reference	4.21		

Only slight amounts of bias between the TSH AccuBind® ELISA method and the reference method are indicated by the closeness of the mean values. The least square regression equation and correlation coefficient indicates excellent method agreement.

##### 14.4 Specificity

The cross-reactivity of the TSH AccuBind® ELISA test system to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of thyrotropin needed to produce the same absorbance.

Substance	Cross Reactivity	Concentration
Thyrotropin (hTSH)	1.0000	-
Foliotropin (hFSH)	< 0.0001	1000ng/ml
Lutropin Hormone (hLH)	< 0.0001	1000ng/ml
Chorionic Gonadotropin (hCG)	< 0.0001	1000ng/ml

#### 14.5 Correlation between 1 hr and 2 hr incubation

The one- (1) hr and two (2) hr (optional) incubation procedures were compared. Thirty (30) biological specimens (ranging from 0.1 - 18.5  $\mu\text{IU/ml}$ ) were used. The least square regression equation and the correlation coefficient were computed for the 2 hr procedure (y) in comparison with the 1 hr method (x). Excellent agreement is evidenced by the correlation coefficient, slope and intercept:  $Y = 0.986 (x) + 0.119$  Regression Correlation = 0.998

#### 15.0 REFERENCES

- Hopton MR, & Harrap JJ, "Immunoradiometric assay of thyrotropin as a first line thyroid function test in the routine laboratory", *Clinical Chemistry*, 32, 691 (1986).
- Caldwell, G et al, "A new strategy for thyroid function testing", *Lancet*, 1, 1117 (1985).
- Young DS, Pestaner LC, and Gilberman U, "Effects of Drugs on Clinical Laboratory Tests", *Clinical Chemistry*, 21, 3660 (1975).
- Spencer, CA, et al, "Interlaboratory/Intermethod differences in Functional Sensitivity of Immunometric Assays of Thyrotropin (TSH) and Impact on Reliability of Measurement of Subnormal Concentrations of TSH", *Clinical Chemistry*, 41, 367 (1995).
- Beck-Peccoz P, Persani L, "Variable biological activity of thyroid stimulating hormone", *Eur J Endocrinol*, 131, 331-340 (1994).
- Bravermann, LE, "Evaluation of thyroid status in patients with thyrotoxicosis", *Clin Chem*, 42, 174-181 (1996).
- Fisher, DA, "Physiological variations in thyroid hormones. Physiological and pathophysiological considerations", *Clin Chem*, 42, 135-139 (1996).

Revision: 4 Date: 2019-Jul-16 DCO: 1353  
MP325 Product Code: 325-300

Size	96(A)	192(B)	480(D)	960(E)
A)	1ml set	1ml set	2ml set	2ml set x2
B)	1 (13ml)	2 (13ml)	1 (60ml)	2 (60ml)
C)	1 plate	2 plates	5 plates	10 plates
D)	1 (20ml)	1 (20ml)	1 (60ml)	2 (60ml)
E)	1 (7ml)	2 (7ml)	1 (30ml)	2 (30ml)
F)	1 (7ml)	2 (7ml)	1 (30ml)	2 (30ml)
G)	1 (8ml)	2 (8ml)	1 (30ml)	2 (30ml)

For Orders and Inquiries, please contact


**Monobind Inc.**  
100 North Pointe Drive  
Lake Forest, CA 92530 USA

Tel: +1 949.951.2685 Mail: info@monobind.com  
Fax: +1 949.951.3539 Fax: www.monobind.com

Please visit our website to learn more about our products and services.

**Glossary of Symbols**  
(EN 980/ISO 15223)

## Anti-Müllerian Hormone (AMH) Test System

Product Code: 9725-300

### 1.0 INTRODUCTION

Intended Use: The Quantitative Determination of AMH Concentration in Human Serum or Plasma by a Microplate Enzyme Immunoassay, Colorimetric

### 2.0 SUMMARY AND EXPLANATION OF THE TEST

Anti-Müllerian Hormone (AMH) is a disulfide-linked homodimeric 140kDa glycoprotein from the trans-forming growth factor- $\beta$  (TGF- $\beta$ ) superfamily of growth and differentiation factors. It is primarily produced by the gonads in both males and females.

In fetal males, AMH is produced by the Sertoli cells and induces regression of the Müllerian duct and therefore promotes development of the male reproductive tract.<sup>1</sup> Infant males have very high levels of AMH (>30 ng/ml) that slowly decreases until post-pubescence where it remains at a low level (<10 ng/ml).<sup>2,3</sup>

In females, AMH begins to be produced near the time of birth with levels increasing until puberty. After puberty, blood AMH levels decrease until menopause where it becomes nearly undetectable (<0.1 ng/ml). AMH concentration in female blood has repeatedly been linked to ovarian reserve, thereby giving an indication to patients' remaining reproductive lifespans.<sup>4</sup> Additionally, high levels of AMH (>4.7 ng/ml, 80% CI) in females are an indication of polycystic ovarian syndrome (PCOS).<sup>4</sup>

When AMH levels drop below 1.0 ng/ml in females, they are considered to have low ovarian reserves. Patients in these ranges are advised to not delay family planning or to undergo infertility treatments such as in vitro fertilization.<sup>4,5</sup>


The AccuBind<sup>®</sup> AMH test kit is a highly sensitive assay that can be used to measure blood AMH levels in order to monitor progress of patients' infertility treatments and approximate the onset of menopause.

### 3.0 PRINCIPLE

#### Sandwich Equilibrium Method (TYPE 2):

The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitope recognition, in excess, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of x-AMH antibody coated on the well.

Upon mixing the enzyme-labeled x-AMH antibody (separate epitope) and serum containing the native antigen, a reaction results between the native antigen and the antibodies without competition or steric hindrance to form a sandwich complex. The interaction is illustrated by the following equation:

After sufficient time results, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

### 4.0 REAGENTS

#### Materials Provided:

**A. AMH Calibrators – 1.0 ml/vial (Lyophilized) Icons [A – F]**  
 Six (6) vials of references for AMH at levels of 0(A), 0.2(B), 0.5(C), 1.0(D), 5.0(E) and 15.0(F) ng/ml. Store at 2-8 °C. **Reconstitute each vial with 1.0ml of distilled or deionized water.** The reconstituted calibrators are stable for 10 days at 2-8 °C. To store for a longer period, aliquot the reconstituted calibrators into cryo vials and store at -20 °C. **DO NOT FREEZE/ THAW MORE THAN TWICE.** A preservative has been added.  
**Note:** The calibrators, human serum based, were calibrated using a reference preparation and are traceable against NIBSC code 16/190.

**B. AMH Controls – 1.0 ml/vial (Lyophilized) Icons [M&N]**  
 Two (2) vials of reference controls for AMH. Store at 2-8 °C. **Reconstitute each vial with 1.0ml of distilled or deionized water.** The reconstituted controls are stable for 10 days at 2-8 °C. To store for a longer period, aliquot the reconstituted calibrators into cryo vials and store at -20°C. **DO NOT FREEZE/ THAW MORE THAN TWICE.** A preservative has been added.

**C. AMH Enzyme Reagent – 6 ml/vial – Icon**   
 One (1) vial contains anti-AMH conjugate reagent. Store at 2-8 °C.

**D. AMH Antibody Coated Plate – 96 wells – Icon**   
 One 96-well microplate coated with x-AMH antibody. Store at 2-8 °C.

**E. Wash Solution Concentrate – 20 ml/vial – Icon**   
 One (1) vial contains a surfactant in buffered saline. A preservative has been added. Store at 2-8 °C. See *Reagent Preparation section*.

**F. Substrate Reagent – 12 ml/vial – Icon**   
 One (1) vial contains tetramethylbenzidine (TMB) and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) in buffer. Store at 2-8 °C.

**G. Stop Solution – 8 ml/vial – Icon**   
 One (1) vial contains a strong acid (0.5 M H<sub>2</sub>SO<sub>4</sub>). Store at 2-8°C.

#### H. Product Instructions.

**Note 1:** Do not use reagents beyond the kit expiration date.

**Note 2:** Do not expose reagents to heat, sun, or strong light.

**Note 3:** The above components are for one 96-well microplate. For other kit configurations, refer to table at the end of the instructions.

#### 4.1 Required But Not Provided:

- Pipette capable of delivering 0.050ml (50µl) and 0.100ml (100µl) volumes with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.100ml (100µl) and 0.350ml (350µl) volumes with a precision of better than 1.5%.
- Microplate washers or a squeeze bottle (optional).
- Microplate reader with 450nm and 620nm wavelength absorbance capability.
- Absorbent paper for blotting the microplate wells.
- Plastic wrap or microplate cover for incubation steps.
- Vacuum aspirator (optional) for wash steps.
- Timer.
- Quality control materials.

### 5.0 PRECAUTIONS

**For In Vitro Diagnostic Use**  
**Not for Internal or External Use in Humans or Animals**

All products that contain human serum have been found to be non-reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA licensed reagents. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-8395.

**Safe Disposal of kit components must be according to local regulatory and statutory requirement.**

### 6.0 SPECIMEN COLLECTION AND PREPARATION

The specimens shall be blood serum or plasma in type and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain redtop venipuncture tube without additives or anti-coagulants for serum or EDTA/heparin containing tubes for plasma. Allow the blood to clot for serum samples. Centrifuge the specimen to separate the serum or plasma from the cells.

If the specimen(s) cannot be assayed immediately after blood withdrawal, the sample(s) may be stored at temperatures of 2-8 °C for up to seven (7) days or -20 °C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing (a maximum of two freeze/thaw cycles prior to use). When assayed in duplicate, 0.100 ml (100 µl) of the specimen is required.

### 7.0 QUALITY CONTROL

Each laboratory should assay controls at levels in the low, medium and high ranges of the dose response curve for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

### 8.0 REAGENT PREPARATION

#### 1. Wash Buffer

Dilute contents of wash solution to 1000ml with distilled or deionized water in a suitable storage container. Diluted buffer can be stored at 2-30°C for up to 60 days.

**Note: Do not use reagents that are contaminated or have bacterial growth.**

### 9.0 TEST PROCEDURE

*Before proceeding with the assay, bring all reagents, serum reference calibrators and controls to room temperature (20-27°C).*

**\*\*Test Procedure should be performed by a skilled individual or trained professional\*\***

- Format the microplates' wells for each serum reference, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8 °C.**
- Pipette 0.050 ml (50 µl) of the appropriate serum reference calibrator, control or specimen into the assigned well.
- Add 0.050 ml (50 µl) of the AMH Enzyme Reagent to each well. **It is very important to dispense all reagents close to the bottom of the coated well.**

- Swirl the microplate gently for 20-30 seconds, cover and incubate for 60 minutes at room temperature.
- Discard the contents of the microplate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
- Add 0.350 ml (350 µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100 µl) of TMB Substrate to all wells. **Always add reagents in the same order to minimize reaction time differences between wells.**  
**DO NOT SHAKE PLATE AFTER SUBSTRATE ADDITION**
- Incubate at room temperature for twenty (20) minutes.
- Add 0.050 ml (50 µl) of stop solution to each well and mix gently for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
- Read the absorbance in each well at 450nm (using a reference wavelength of 630nm to minimize well imperfections) in a microplate reader. **The results should be read within fifteen (15) minutes of adding the stop solution.**

Note: For re-assaying specimens with concentrations greater than 15 ng/ml, the samples should be diluted and multiplied accordingly. Acceptable diluents are serum from postmenopausal females (<0.1 ng/ml AMH), the "0" calibrator, and other diluent solutions sold by Monobind.

### 10.0 CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of AMH in unknown specimens.


- Plot the absorbance for each duplicate serum reference versus the corresponding AMH concentration in ng/ml on linear graph paper (do not average the duplicates of the serum references before plotting).
- Draw the best-fit curve through the plotted points.
- To determine the concentration of AMH for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (0.299) intersects the dose response curve at 1.33 ng/ml AMH concentration (See Figure 1).

**Note:** Computer data reduction software designed for ELISA assays may also be used for the data reduction. **If such software is utilized, the validation of the software should be ascertained.**

EXAMPLE 1

Sample I.D.	Well Number	Abs (A)	Mean Abs (B)	Value (ng/ml)
Cal A	A1	0.005	0.005	0
	B1	0.005		
Cal B	C1	0.041	0.043	0.20
	D1	0.045		
Cal C	E1	0.109	0.112	0.50
	F1	0.115		
Cal D	G1	0.215	0.225	1.00
	H1	0.234		
Cal E	A2	1.047	1.068	5.00
	B2	1.097		
Cal F	C2	2.570	2.635	15.00
	D2	2.678		
Patient	E2	0.300	0.299	1.33
	F2	0.298		

\*The data presented in Example 1 and Figure 1 is for illustration only and **should not** be used in lieu of a dose response curve prepared with each assay.

\*If the absorbance readout is off-scale or higher than the average absorbance of the highest calibrator, sample should be repeated with dilution.

### 11.0 Q.C. PARAMETERS

In order for the assay results to be considered valid the following criteria should be met:

1. Maximum Absorbance (Calibrator 'F')  $\geq 1.3$
2. Four out of six quality control pools should be within the established ranges.

### 12.0 RISK ANALYSIS

The MSDS and Risk Analysis Form for this product are available on request from Monobind Inc.

#### 12.1 Assay Performance

1. It is important that the time of reaction in each well is held constant to achieve reproducible results.
2. Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
3. Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
4. If more than one (1) plate is used, it is recommended to repeat the dose response curve.
5. The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
6. Plate readers measure vertically. Do not touch the bottom of the wells.
7. Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
8. Use components from the same lot. No intermixing of reagents from different batches.
9. Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind's IFU may yield inaccurate results.
10. All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
11. It is important to calibrate all the equipment e.g. pipettes, readers, washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
12. Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from [Monobind@monobind.com](mailto:Monobind@monobind.com).

#### 12.2 Interpretation

1. **Measurements and interpretation of results must be performed by a skilled individual or trained professional.**
2. Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
3. For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
4. If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability.**
5. If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

6. The reagents for AccuBind® ELISA procedure have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays. (Boscato LM Stuart MC. *Heterophilic antibodies: a problem for all immunoassays* Clin. Chem 1988;34:27-33). For diagnostic purposes the results from this assay should be used in combination with clinical examination, patient's history, and, all other clinical findings.

### 13.0 EXPECTED RANGES OF VALUES

AMH levels were measured by the AMH AccuBind® Test System in apparently normal females of different age groups. The values obtained are given in Table 2.

**Table 2**  
Female Reference Ranges for the AMH Test System

Age Group	N	Mean (ng/ml)	Median (ng/ml)	Minimum (ng/ml)	Maximum (ng/ml)
20-29	24	5.95	5.94	1.77	12.41
30-39	80	2.83	2.47	0.11	12.67
40-49	38	1.33	0.59	0.02	9.77

The expected values for males and females of untested age groups as determined by literature sources are collected in Table 3.

**Table 3**  
Additional Reference Ranges for the AMH Test System

Males	Range (ng/ml)
<24 months	14-466
2-12 years	7.4-243
>12 years	0.7-19
Females	
<24 months	<4.7
2-12 years	<8.8
13-19 years	0.9-9.5

It is important to keep in mind that establishment of a range of values, which can be expected to be found by a given method for a population of "normal" persons, is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons, each laboratory should depend upon the range of expected values established by the manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

### 14.0 PERFORMANCE CHARACTERISTICS

#### 14.1 Precision

The within assay precision of the AMH AccuBind® Microplate EIA Test System were determined by analyses on six different levels of pool control and patient sera. The mean values, standard deviation and coefficient of variation for each of these control sera are presented in Table 4.

**TABLE 4**  
Precision data for the AMH Test System

Sample	Mean Value (pg/ml)	Within-Run Precision		Total Precision (n=80)	
		SD	CV%	SD	CV%
Control 1	0.80	0.02	2.87	0.04	5.32
Control 2	4.54	0.12	2.60	0.24	5.24
Control 3	14.16	0.30	2.11	0.77	5.47
Patient 1	1.42	0.03	2.20	0.08	5.79
Patient 2	0.30	0.02	5.69	0.03	9.46
Patient 3	9.57	0.26	2.73	0.51	5.36

\*As measured in forty experiments in duplicate over a 20 day period.


#### 14.2 Sensitivity

The AMH AccuBind® ELISA test system has a LoB=0.024 ng/ml and a LoD=LoQ=0.044 ng/ml.

#### 14.3 Accuracy

#### 14.3.1 Linearity

The linearity of the AMH AccuBind® ELISA test system was tested by diluting human serum samples containing high levels of AMH (2 to 17.5 ng/ml) with low AMH (<0.2 ng/ml) human serum samples or the "0" Calibrator. The system demonstrates excellent linearity through the range of the test up to 17.5 ng/ml as shown in Figure 2.


#### 14.3.2 Recovery

The recovery of the AMH AccuBind® Microplate ELISA Test System was calculated for five patient samples spiked with 0.5, 1.0, 2.0, 4.0, and 8.0 ng/ml AMH. Recoveries were determined to be within 15% of the expected values for all samples.

#### 14.3.3 Method Comparison

The AMH AccuBind® Microplate ELISA Test System was initially evaluated on 167 patients with known AMH concentrations calculated by a different AMH test system from a different manufacturer. Correlation between the two methods is excellent with a R<sup>2</sup> coefficient = 0.951. A graph of the data is shown in Figure 3.

#### 14.4 Cross-Reactivity

The following analytes were tested and found to be non-reactive.

Analyte	Concentration	% Reactivity
Follicle Stimulating Hormone	100 mIU/ml	<0.001
Human Chorionic Gonadotropin	1000 mIU/ml	<0.001
Leuteinizing Hormone	200 mIU/ml	<0.001
Prolactin	100 ng/ml	<0.001

### 15.0 REFERENCES

1. Wilson, J.D.; George, F.W.; Griffin, J.E. The hormonal control of sexual development. *Science* 1981, 211, 4488, 1278-1284.
2. Jopling, H.; Yates, A.; Burgoyne, N.; Hayden, K.; Chaloner, C.; Tetlow, L. Paediatric Anti-Müllerian Hormone measurement: Male and Female reference intervals established... *Endocrinol Diab Metab*. 2018, 1, 4, e00021.
3. Matuszczak, E.; Hermanowicz, A.; Komarowska, M.; Debek, W. Serum AMH in Physiology and Pathology of Male Gonads. *Int. J. of Endocrin.* 2013, 2013, Article ID 128907, 6 pgs.
4. Broer, S.L.; Broekmans, F.J.M.; Laven, J.S.E.; Fauser, B.C.J.M. Anti-Müllerian hormone: ovarian reserve testing and its potential clinical implications. *Human Reproduction Update*, 2014, 20, 5, 688-701.
5. Lee, J.Y.; Ahn, S.; Lee, J.R.; Jee, B.C.; Kim, C.H.; Seo, S.; Suh, C.S.; Kim, S.H. Reference Values for the Revised Anti-Müllerian Hormone Generation II Assay: Infertile Population-based Study. *J Korean Med Sci* 2017, 32, 825-829.
6. Test ID: AMH Antimüllerian Hormone (AMH), Serum. <https://www.mayocliniclabs.com/test-catalog/Clinical+and+Interpretive/89711> (accessed Aug 16, 2019).

Size	96(A)	192(B)	
Reagent (fill)	A)	1.0ml (dry) set	1.0ml (dry) set
	B)	1.0ml (dry) set	1.0ml (dry) set
	C)	1 (6ml)	2 (6ml)
	D)	1 plate	2 plates
	E)	1 (20ml)	1 (20ml)
	F)	1 (12ml)	2 (12ml)
	G)	1 (8ml)	2 (8ml)

For Orders and Inquires, please contact


**Monobind Inc.**  
100 North Pointe Drive  
Lake Forest, CA 92630 USA

Tel: +1 949.951.2665 Mail: [info@monobind.com](mailto:info@monobind.com)  
Fax: +1 949.951.3539 Fax: [www.monobind.com](http://www.monobind.com)

Please visit our website to learn more about our products and services.

### Glossary of Symbols (EN 980/ISO 15223)

## Instructions for Use

# Sperm Antibody ELISA

IVD

CE

REF EIA-1826

 96

**DRG Instruments GmbH**, Germany  
Frauenbergstraße 18, 35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

Distributed by:

**DRG International, Inc.**, USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.**  
**Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.**  
**Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.**  
**Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.**  
**Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**  
**Utilize apenas a versão válida das Instruções de Utilização fornecidas com o kit.**

<b>Introduced modifications / Durchgeführte Änderungen / Modifiche introdotte / Modificaciones introducidas / Modifications apportées / Modificações introduzidas</b>	
<p>The following changes have been made in comparison to the previous version:            Im Vergleich zur Vorgängerversion wurden folgende Änderungen vorgenommen:            Rispetto alla versione precedente, sono state apportate le seguenti modifiche:            Se han introducido los siguientes cambios en comparación con la versión anterior:            Les modifications suivantes ont été apportées par rapport à la version précédente :            Foram efetuadas as seguintes alterações em comparação com a versão anterior:</p>	
<p><b>Detailed editorial revision.</b> Changed wording in several chapters.  <b>Ausführliche redaktionelle Überarbeitung.</b> Geänderter Wortlaut in mehreren Kapiteln.  <b>Revisione editoriale dettagliata.</b> Modificato il testo in diversi capitoli.  <b>Revisión editorial detallada.</b> Se ha cambiado la redacción de algunos capítulos.  <b>Révision éditoriale détaillée.</b> Modification de la formulation dans plusieurs chapitres.  <b>Revisão editorial detalhada.</b> Texto alterado em vários capítulos.</p>	
1 INTENDED USE:	More definitions for intended use and use not intended.
4.4 Reagent Preparation:	Correction for stability of diluted wash solution to 1 week (before: 2 weeks).
5.2 Samples Storage:	Sample stability at 2 °C to 8 °C changed to 2 months (before: 12 months)
6.2 Test Procedure:	<p><b>Incubation time for Enzyme Conjugate</b> reduced to <b>30 minutes</b>            (before: 60 minutes);  <b>Pipetting volume for Stop Solution</b> increased to <b>100 µL</b>            (before: 50 µL);</p>
9 PERFORMANCE CHARACTERISTICS:	Updated and additional data.
10 LIMITATIONS OF THE PROCEDURE:	Updated.

ENGLISH .....	2
DEUTSCH .....	10
ITALIANO .....	17
ESPAÑOL .....	24
FRANÇAIS .....	31
LITERATURE .....	38
SYMBOLS USED .....	39

## 1 INTENDED USE

The **DRG Sperm Antibody ELISA** is a manual enzyme immunoassay for the **quantitative** measurement of antibodies directed against human spermatozoa (ASA) in human serum.

**For *in vitro* diagnostic use. For laboratory professional use.**

The device is **intended** to be used as an aid for the diagnosis of immunologically caused disorders of fertility.

The device is **not intended** to be used for the detection of poly- and monoclonal gammopathies.

### 1.1 Scientific Validity

Antibodies directed against spermatozoa antigens may cause infertility in women or men. The application of the Sperm Antibody ELISA is recommended for the diagnosis of immunologically caused disorders of fertility.

Unwanted childlessness is a growing problem with which up to 20% of all couples in the reproductive age are confronted temporarily or long-term. In 5-20 % of these cases the presence of anti-spermatozoa antibodies in the male or the female patient is detectable [1,2,15].

The definition of infertility according to the WHO (WHO Laboratory Manual for the Examination of Human Semen and Semen Cervical-Mucus Interaction, 1999) is the absence of a conception within 12 months of unprotected intercourse. The main cause of an immunological fertility disorder is the formation of antibodies directed against spermatozoa antigens.

Anti-spermatozoa antibodies (ASA) exert heterogeneous effects on the ability of spermatozoa to fertilize. The inhibiting effect of ASA on the motility of spermatozoa by binding to their surface and by agglutinating processes is well-known [3].

The penetration of the spermatozoa into the cervical mucus is impaired by the presence of ASA in the seminal plasma and/or in the cervical mucus [4]. ASA negatively influence the capacitation and the acrosome reaction of spermatozoa and thereby impede the interaction of the spermatozoa with the oocyte [5,6].

The interaction of the spermatozoon with the oocyte and the subsequent binding to and penetration of the zona pellucida may be inhibited by ASA. The following fusion of the oocyte and a spermatozoon may also be impaired by the presence of ASA [7,8].

The rate of pregnancies in couples with ASA on the part of the man or the woman was shown to be 38% lower compared to the control groups [9]. Furthermore, an influence on the implantation and on the early embryological development could be confirmed. An association of ASA and miscarriages is discussed. The frequency of ASA in infertile couples amounts to 20% [10,11].

ASA may occur dissolved in the ejaculate or bound to the surface of spermatozoa. In addition, ASA can be detected in serum. ASA may be found in men and in women [12]. In women, ASA may be found in cervical mucus, oviduct liquid and follicular liquid. Men having more than 50% of their spermatozoa coated with anti-spermatozoa antibodies show a conspicuously reduced rate of fertility [13].

ASA have been shown to be associated with chronic prostatitis which has a negative effect on male reproductive function [14].

## 2 PRINCIPLE OF THE TEST

The DRG Sperm Antibody ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **sandwich principle**.

The microtiter wells are coated with a mix of spermatozoa proteins.

During incubation, anti-spermatozoa antibodies in the samples (standards, controls, patient samples) bind to the coated surface of the wells.

A washing step removes unbound sample components.

Added enzyme conjugate binds to the immobilized antigen-antibody-complexes. The conjugate contains anti-human immunoglobulin antibodies, labelled with horseradish peroxidase (HRP).

After a washing step to remove all unbound substances, the solid phase is incubated with the substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured. The intensity of color is proportional to the concentration of the analyte in the sample.

A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

### 3 WARNINGS AND PRECAUTIONS

- This kit is for *in vitro* diagnostic use only. For laboratory professional use only.
- Before starting the assay, read the instructions for use completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
- Do not mix or use components from kits with different lot numbers. It is advised not to interchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
- Do not use reagents beyond expiry date as shown on the kit labels.
- Do not reuse microtiter wells.
- Reagents of other manufacturers must not be used together with the reagents of this test kit.
- All reagents in this kit are clear liquids, substrate solution is clear and colorless. Changes in its appearance may affect the performance of the test. In that case, contact DRG.
- Microbial contamination of reagents or samples may give false results.
- Allow the reagents to reach room temperature (20 °C to 26 °C) before starting the test. Temperature will affect the optical density readings of the assay.
- All indicated volumes must be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
- Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into original vials as reagent contamination may occur.

#### General precautions

- Follow good laboratory practice and safety guidelines.
- Never pipet by mouth and avoid contact of reagents and samples with skin and mucous membranes.
- Do not smoke, eat, drink, or apply cosmetics in areas where samples or kit reagents are handled.
- Wear lab coats and disposable latex gloves when handling samples and reagents and where necessary safety glasses.

#### Biohazard information

- All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. However, no known test method can offer total assurance that no infectious agent is present.
- The device contains material of animal origin, which is certified apparently free of infectious or contagious diseases and injurious parasites.
- Bovine components originate from countries where BSE (Bovine spongiform encephalopathy) has not been reported.
- All materials and samples of human or animal origin must be handled as if capable of transmitting infectious diseases.
- Handling must be done in accordance with the procedures defined by appropriate national biohazard and safety guideline or regulation. Waste must be discarded according to local rules and regulations.

#### Information to chemical hazards and hazard classification

- Some reagents contain preservatives in non-declarable concentrations. Nevertheless, in case of contact with eyes or skin, flush immediately with water.
- Substrate Solution contains an ingredient in non-declarable concentrations which causes serious eye irritation. In case of possible contact with eyes, rinse immediately carefully and thoroughly with eye wash or water. After contact with skin, wash with plenty of water. Take-off contaminated clothing and wash it before reuse.
- Avoid contact with Stop Solution containing 0.5 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.
- Chemicals and prepared or used reagents must be treated as hazardous waste according to the national safety guideline or regulation.
- This product does not contain substances which have carcinogenic, mutagenic or toxic for reproduction (CMR) properties.

All reagents of this test kit do NOT contain hazardous substances in concentrations to be declared, a classification and labelling is not required.

For detailed information, please refer to the Safety Data Sheet, which is available upon request directly from DRG.

## 4 MATERIALS

### 4.1 Materials provided with the kit

Symbol	Quantity	Description	Preparation
<b>Microtiterwells</b>	12 x 8 wells (break apart)	<b>Microtiter plate</b> Coated with a mix of spermatozoa proteins	Ready to use
<b>Dilution Buffer / Zero Standard</b>	1 x 50 mL	<b>Sample Dilution Buffer / Zero Standard *</b> Concentration: 0 U/mL	Ready to use
<b>Standard (Standard 1 - 4)</b>	4 x 0.5 mL	<b>Standards *</b> Concentrations: 31 – 62 – 125 – 250 U/mL	Ready to use
<b>Quality Control</b>	1 x 0.5 mL	<b>Control *</b> <i>For control values and ranges please refer to vial label or Certificate of Analysis.</i>	Ready to use
<b>Enzyme Conjugate</b>	1 x 8 mL	<b>Enzyme Conjugate *</b> Anti-human IgG antibody conjugated to horseradish peroxidase; Colored red.	Ready to use
<b>Substrate Solution</b>	1 x 14 mL	<b>Substrate Solution</b> Contains 3,3',5,5'-tetramethylbenzidine (TMB). <i>Keep away from direct sun light.</i>	Ready to use
<b>Stop Solution</b>	1 x 14 mL	<b>Stop Solution</b> Contains 0.5 M H <sub>2</sub> SO <sub>4</sub> . <i>Avoid contact with the stop solution. It may cause skin irritations and burns.</i>	Ready to use
<b>Wash Solution</b>	1 x 30 mL	<b>Wash Solution, 40X concentrate *</b>	See "Reagent Preparation".
<b>Cover foil</b>	1 x	<b>Cover foil</b>	
	1 x	<b>Instructions for Use</b>	
	1 x	<b>Certificate of Analysis (CoA)</b>	
		* Contains non-mercury preservative.	

### 4.2 Materials required but not provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm)
- Calibrated variable precision micropipettes
- Incubator for 37 °C
- Manual or automatic equipment for rinsing microtiter plate wells
- Absorbent paper
- Distilled water
- Timer
- Graph paper or software for data reduction

### 4.3 Storage and Stability of the Kit

**Unopened kits and reagents** as well as **opened reagents** must be stored at 2 °C to 8 °C.

The microtiter plate contains snap-off strips. Do not open the pouch of the wells until it reaches room temperature. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch including the desiccant and used in the plate frame provided. Once the foil bag has been opened, care must be taken to close it tightly again.

Once opened, reagent vials must be closed tightly again.

	Storage Temperature	Stability
Unopened kits and unopened reagents	2 °C to 8 °C	Until the expiration date printed on the label. Do not use reagents beyond this date!
Opened kit	2 °C to 8 °C	8 weeks

#### 4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (20 °C to 26 °C) prior to use.

##### Wash Solution

Add distilled water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL distilled water to a final volume of 1200 mL.

Stability after dilution:	at 20 °C to 26 °C	1 week
---------------------------	-------------------	--------

#### 4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

#### 4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components must not be used for a test run. They have to be stored until a final solution has been found. After this, they must be disposed of according to the official regulations.

### 5 SAMPLE COLLECTION, STORAGE AND PREPARATION

The following sample material can be used in this test:

#### Human serum

Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use hemolytic, icteric, or lipemic samples. For further information refer to chapter “*Interfering Substances*”.

*For the determination of anti-spermatozoa antibodies in seminal plasma please use our Sperm Antibody (seminal plasma) ELISA (REF EIA-4249).*

#### 5.1 Sample Collection

**Serum:** Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Whole blood should not be frozen before centrifugation.

#### 5.2 Samples Storage

Samples must be stored tightly capped prior to performing the assay. If stored frozen, freeze only once. Thawed samples must be inverted several times prior to testing.

Stability	at 2 °C to 8 °C	4 days
	at -20 °C (in aliquots)	up to 2 months

#### 5.3 Sample Preparation

Prior to assaying, dilute each patient sample **1:100** with *Dilution Buffer*.

##### Example:

Dilution 1:100: 5 µL sample + 495 µL *Dilution Buffer* (mix thoroughly)

Stability of diluted samples	at 2 °C to 8 °C	4 days
	at -20 °C (in aliquots)	7 days

**Note:** The *Quality Control* is ready to use and must not be diluted!

### 6 ASSAY PROCEDURE

#### 6.1 Procedural Notes

- All reagents and samples must be allowed to come to room temperature (20 °C to 26 °C) before use.
- All reagents must be mixed without foaming.
- Do not interchange caps of reagent vials to avoid cross-contamination.
- Use new disposal plastic pipette tips for each standard, control, or sample in order to avoid carry-over.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense conjugate without splashing accurately to the bottom of wells.
- Mix the contents of the microtiter plate wells thoroughly to ensure good test results.

- Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
- Once the test has been started, all steps must be completed without interruption and in the same sequence for each step.
- The enzymatic reaction is linearly proportional to time and temperature.
- Optical density is a function of the incubation time and temperature. Respect the incubations times and temperatures as given in chapter "Test Procedure".
- Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- During the incubation at 37 °C cover microtiter strips with foil to avoid evaporation.
- **Important note to wash procedure:**  
Washing is critical. Improperly washed wells will give erroneous results. The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
- **Test performance using fully automated analysis devices:**  
Automated test performance using fully automated, open-system analysis devices is possible. However, the combination must be validated by the user.

## 6.2 Test Procedure

Each run must include a standard curve.

The controls serve as internal controls for the reliability of the test procedure. They must be assayed with each test run.

The given test procedure describes manual processing.

**Important note:** The accuracy of this assay is markedly influenced by the correct incubation temperature and incubation time. Pipetting the samples should not exceed 15 minutes.

1. Secure the desired number of microtiter wells in the frame holder.
2. Pipette **50 µL** of each **Zero Standard, Standard, Quality Control**, and **diluted sample with new disposable tips** into appropriate wells.
3. Cover with foil and incubate for **60 minutes at 37 °C**.
4. Wash the wells as follows:  
If the wash step is performed manually:  
Briskly shake out the contents of the wells.  
Rinse the wells **3 times** with **300 µL** diluted *Wash Solution* per well.  
  
If an automated plate washer is used:  
Rinse the wells **3 times** with **400 µL** diluted *Wash Solution* per well.  
  
At the end of the washing step, always strike the wells sharply on absorbent paper to remove residual droplets!
5. Add **50 µL Enzyme Conjugate** into each well.
6. Cover with foil and incubate for **30 minutes at 37 °C**.
7. Wash as described in step 4.
8. Pipette **50 µL of Substrate Solution** to each well.
9. Incubate for **30 minutes** at room temperature.
10. Stop the enzymatic reaction by adding **100 µL of Stop Solution** to each well.
11. Measure the optical density (OD) of the solution in each well at **450 nm (reading) and at 620 nm to 630 nm (background subtraction, recommended)** with a microtiter plate reader.  
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

## 6.3 Calculation of Results

1. The concentration of the samples can be read directly from the standard curve.  
The standards are already pre-diluted, therefore the 1:100 dilution of the samples must not be taken into account for the final calculation of sample concentrations.
2. For duplicate determinations, the mean of the two optical density (OD) values for each standard, control, and patient sample must be taken. If the two values deviate substantially from one another, DRG recommends retesting the samples.
3. Samples with concentrations exceeding the highest standard can be further diluted with *Dilution Buffer* and re-assayed as described in "Test Procedure", or must be reported as > 250 U/mL. For the calculation of the concentrations, this dilution factor must be considered.
4. Automated method:  
The results in the instructions for use have been calculated automatically using a four-parameter logistic (4PL) curve fit. (4PL Rodbard or 4PL Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. Manual method:  
Using graph paper, construct a standard curve by plotting the (mean) OD obtained from each standard against its concentration with OD value on the vertical (Y) axis and concentration on the horizontal (X) axis.  
Determine the corresponding sample concentration from the standard curve by using the (mean) OD value for each sample.

### 6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Density (450 nm)
<i>Zero Standard</i> (0 U/mL)	0.147
<i>Standard 1</i> (31 U/mL)	0.515
<i>Standard 2</i> (62 U/mL)	0.857
<i>Standard 3</i> (125 U/mL)	1.423
<i>Standard 4</i> (250 U/mL)	2.127

## 7 REFERENCE VALUES

It is strongly recommended that each laboratory determine its own reference values.

In a study conducted with apparently healthy subjects, using the DRG Sperm Antibody ELISA the following data were observed:

Normal values	0 U/mL – 60 U/mL
Borderline	55 U/mL – 65 U/mL
Elevated values	> 60 U/mL

In the case of a value in the range near the cut-off (55 U/mL to 65 U/mL) we recommend a follow-up determination using a new sample taken within the next two weeks.

The results alone should not be the only reason for any therapeutic consequences. The results must be correlated to other clinical observations and diagnostic tests.

## 8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day-to-day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the Quality Control Laboratory are stated in the Certificate of Analyses (CoA) added to the kit. The values and ranges stated on the CoA always refer to the current kit lot and must be used for direct comparison of the results.

If available, it is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Apply appropriate statistical methods for analyzing control values and trends. If the results of the assay do not agree with the established acceptable ranges of control materials, patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above-mentioned items without finding any error contact your distributor or DRG directly.

## 9 PERFORMANCE CHARACTERISTICS

### 9.1 Sensitivity

Limit of Blank (LoB)	0.490 U/mL
Limit of Detection (LoD)	3.367 U/mL
Limit of Quantification (LoQ)	9.632 U/mL
Measuring range	3.367 U/mL – 250 U/mL
Linear range	4.333 U/mL – 250 U/mL

### 9.2 Reproducibility

#### 9.2.1 Within-run Precision

The within-run precision was determined with 4 patient samples covering the complete measuring range within 20 days in 2 independent runs per day. CV was calculated as mean CV of 40 runs.

Sample	n	Mean (U/mL)	CV (%)
1	40	39.68	3.9
2	40	77.99	3.0
3	40	106.16	4.2
4	40	155.95	3.9

**9.2.2 Between-run Precision**

The between-run precision was determined with 4 patient samples covering the complete measuring range within 20 days in 2 independent runs per day and with 2 replicates per run (20 x 2 x 2). CV was calculated from 80 determinations.

Sample	n	Mean (U/mL)	CV (%)
1	80	39.68	11.2
2	80	77.99	8.7
3	80	106.16	11.4
4	80	151.86	13.1

**9.2.3 Between-lot Precision**

The between-lot variation was determined by 6 measurements of different samples with 3 different kit lots.

Sample	n	Mean (U/mL)	CV (%)
1	18	17.78	9.3
2	18	26.65	10.4
3	18	34.58	3.9
4	18	55.56	14.7

**9.3 Recovery**

Recovery was determined by adding increasing amounts of the analyte to different patient samples containing different amounts of endogenous analyte. The percentage recoveries were determined by comparing expected and measured values of the samples.

		Sample 1	Sample 2	Sample 3	Sample 4
<b>Concentration (U/mL)</b>		10.86	34.93	38.34	45.42
<b>Average Recovery (%)</b>		98.6	89.0	101.0	105.0
<b>Range of Recovery (%)</b>	<b>from</b>	90.3	86.2	92.0	99.0
	<b>to</b>	114.4	93.1	111.3	109.5

**9.4 Linearity**

Samples containing different amounts of analyte were serially diluted with *Dilution Buffer*. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

		Sample 1	Sample 2	Sample 3	Sample 4
<b>Concentration (U/mL)</b>		36.37	66.20	183.74	280.32
<b>Average Recovery (%)</b>		105.7	95.1	103.2	104.7
<b>Range of Recovery (%)</b>	<b>from</b>	103.9	92.3	99.7	98.6
	<b>to</b>	107.5	97.9	110.8	107.3

**10 LIMITATIONS OF THE PROCEDURE**

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the instructions for use and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

**10.1 Interfering Substances**

Hemoglobin (up to 4 mg/mL), bilirubin (up to 0.5 mg/mL) and triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

Sera from patients with liver diseases should not be used.

Results may be adversely affected by certain pathologic conditions, such as poly- and monoclonal gammopathies, autoimmune diseases or by an altered immune status.

**10.2 High-Dose Hook Effect**

"High-Dose Hook Effect" is not detected up to 5000 U/mL of ASA.

## **11 LEGAL ASPECTS**

### **11.1 Reliability of Results**

The test must be performed exactly as per the manufacturer's instructions for use. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. If there is any doubt or concern regarding a result, please contact DRG.

### **11.2 Therapeutic Consequences**

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### **11.3 Liability**

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid. Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

### **11.4 Reporting of Serious Incident**

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

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### 11.4 Notification des incidents graves

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## 1 INTENDED USE

The **DRG  $\beta$ -HCG ELISA Kit** is an enzyme immunoassay for the quantitative *in vitro* diagnostic measurement of total human chorionic gonadotropin (hCG and  $\beta$ -hCG) in serum.

**This kit is NOT intended to be used for the risk evaluation of trisomy 21.**

### 1.1 Summary and Explanation

Chorionic Gonadotropin (hCG) is a glycoprotein hormone which is normally produced by the placenta during pregnancy. After conception, the hCG concentration increases rapidly to reach a peak near the end of the first trimester. High concentrations are observed throughout pregnancy. After delivery, hCG levels fall rapidly and become undetectable after a few days.

Structurally intact hCG molecules are composed of an alpha and a beta subunit. The alpha subunit is nearly identical to the alpha subunits of other glycoprotein hormones, such as Thyroid Stimulating Hormone (TSH), Luteinizing Hormone (LH), and Follicle Stimulating Hormone (FSH): The differences in the beta subunit of the respective hormones account for their biological specificity and immunochemical distinctiveness.

Monoclonal antibodies recognizing unique sites on the beta chain of the  $\beta$ -hCG/hCG molecule are essential for differentiation between hCG and LH, FSH and TSH.

Specific assays for beta-hCG permit the early detection of pregnancy.

In addition to the elevated hCG levels during pregnancy, high concentrations of  $\beta$ hCG/hCG may be associated with neoplasms of trophoblastic and nontrophoblastic origin such as hydatiform mole, chorionepithelioma, embryonal cell carcinoma, seminoma and many others.

## 2 PRINCIPLE OF THE TEST

The DRG  $\beta$ -HCG ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA) based on the **sandwich principle**.

**The microtiter wells are coated with a monoclonal antibody [mouse] directed towards a unique antigenic site on a  $\beta$ -HCG molecule.** An aliquot of patient sample containing endogenous  $\beta$ -HCG and/or hCG is incubated in the coated well with enzyme conjugate, which is an anti-  $\beta$ -HCG antibody conjugated with horseradish peroxidase. After incubation the unbound conjugate is washed off.

The amount of bound peroxidase is proportional to the concentration of  $\beta$ -HCG/hCG in the sample.

Having added the substrate solution, the intensity of colour developed is proportional to the concentration of  $\beta$ -HCG/hCG in the patient sample.

### 3 WARNINGS AND PRECAUTIONS

1. This kit is for *in vitro* diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV 1/1, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution colored. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay. However, values for the patient samples will not be affected.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with Stop Solution containing 0.5 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

#### 4 REAGENTS

##### 4.1 Reagents provided

1. *Microtiterwells*, 12 x 8 (break apart) strips, 96 wells;  
Wells coated with anti-β-HCG antibody (monoclonal).
2. *Standard (Standard 0-5)*, 6 vials (lyophilized), 1 mL;  
Concentrations: 0, 5; 25; 50; 100; 200 mIU/mL  
Conversion: 1 pg/mL = 0.00916 mIU/mL  
*The standards are calibrated against the following reference material: 5<sup>th</sup> WHO International Standard Chorionic Gonadotrophin, NIBSC code: 07/364.*  
See „Preparation of Reagents“;  
Contain non-mercury preservative.
3. *Sample Diluent*, 1 vial, 10 mL, ready to use,  
Contains non-mercury preservative.
4. *Enzyme Conjugate*, 1 vial, 11 mL, ready to use,  
Anti-β-HCG antibody conjugated to horseradish peroxidase;  
Contains non-mercury preservative.
5. *Substrate Solution*, 1 vial, 14 mL, ready to use,  
Tetramethylbenzidine (TMB).
6. *Stop Solution*, 1 vial, 14 mL, ready to use,  
contains 0.5 M H<sub>2</sub>SO<sub>4</sub>,  
Avoid contact with the stop solution. It may cause skin irritations and burns.

**Note:** Additional *Sample Diluent* for sample dilution is available upon request.

##### 4.2 Materials required but not provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm)
- Calibrated variable precision micropipettes
- Absorbent paper
- Distilled water
- Timer
- Graph paper or software for data reduction

##### 4.3 Storage Conditions

When stored at 2 °C - 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C - 8 °C. Microtiter wells must be stored at 2 °C - 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for 8 weeks if stored as described above.

##### 4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature prior to use.

##### *Standards*

Reconstitute the lyophilized contents of each vial with 1 ml distilled water and let stand for at least 10 minutes at room temperature. Mix several times before use.

*Note:* The reconstituted standards are stable for 2 months at 2 °C - 8 °C.

*For longer storage aliquot and freeze at -20 °C.*

##### 4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

##### 4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

## 5 SPECIMEN COLLECTION AND PREPARATION

Serum should be used in this assay.

*Note:* Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use haemolytic, icteric, or lipaemic specimens. For further information refer to chapter "Interfering Substances".

### 5.1 Specimen Collection

#### Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

### 5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 7 days at 2 °C - 8 °C prior to assaying.

Specimens held for a longer time (up to 12 months) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

### 5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Sample Diluent* and reassayed as described in *Assay Procedure*.

For the calculation of the concentrations this dilution factor has to be taken into account.

#### Example:

- a) dilution 1:10:            10 µL Serum + 90 µL *Sample Diluent* (mix thoroughly)
- b) dilution 1:100:        10 µL dilution a) 1:10 + 90 µL *Sample Diluent* (mix thoroughly).
- c) dilution 1:1000        10 µL dilution b) 1:100 + 90 µL *Sample Diluent* (mix thoroughly).

#### NOTE:

Sera of pregnant women must be diluted 1/100 in *Sample Diluent* before starting the assay.

## 6 ASSAY PROCEDURE

### 6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Absorbance is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

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1. Secure the desired number of Microtiter wells in the holder.
2. Dispense **25 µL** of each *Standard*, *Control* and **samples** with *new disposable tips* into appropriate wells.
3. Dispense **100 µL Enzyme Conjugate** into each well.  
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Incubate for **60 minutes** at room temperature.
5. Briskly shake out the contents of the wells. Rinse the wells 5 times with distilled water (400 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.  
**Important note:** The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
6. Add **100 µL** of *Substrate Solution* to each well.
7. Incubate for **15 minutes** at room temperature.
8. Stop the enzymatic reaction by adding **50 µL** of *Stop Solution* to each well.
9. Determine the optical density of the solution in each well at **450 nm (reading) and at 620 - 630 nm (background subtraction, recommended)** with a microtiter plate reader.  
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

## 6.3 Calculation of Results (quantitative)

1. Calculate the average optical density (OD) values for each set of standards, controls and patient samples.
2. Using linear graph paper, construct a standard curve by plotting the mean OD obtained from each standard against its concentration with OD value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean OD value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 200 mIU/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

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## 6.4 Assay Procedure (qualitative method)

This procedure differentiates positive (pregnant) from negative samples by comparing the sample beta hCG levels with the *Standard 0* (0 mIU/mL) and *Standard 3* (50 mIU/mL).

Patient samples are run with the *Standard 0* and the 50 mIU/mL *Standard*. The assay procedure is identical with the Quantitative Method, but step 9 and 10 is omitted.

## 6.5 Calculation of Results (qualitative)

For a qualitative analysis of the β-hCG level the color development of the specimen is compared with the color of the 0 mIU/mL and 50 mIU/mL standards.

If the blue color is less intense than the color of the 50 mIU/mL standard, the sample is considered as negative.

If the blue color is more intense than or equal to the color of the 50 mIU/mL standard the sample is considered as positive.

## 7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

### 7.1 Normal healthy adults, non-pregnant

In a study conducted with apparently normal healthy adults, using the DRG  $\beta$ -HCG ELISA the following values are observed:

Population	Age (years)	Valid N	$\beta$ -hCG [mIU/mL]
men	< 50	40	< 5
	> 50	10	< 5
women	< 50	42	< 5
	> 50	7	< 5

### 7.2 Pregnant women in the 2<sup>nd</sup> trimester

Week of pregnancy	Valid samples	5% Percentile [mIU/mL]	95% Percentile [mIU/mL]
14	103	10 303	71 980
15	96	9 246	51 666
16	91	5 266	36 947
17	14	4 632	24 033

#### CAUTION:

- For the detection of pregnancy in serum, a qualitative assay is used with a cut-off point of 50 mIU/mL. Negative or borderline results should be repeated on a fresh specimen obtained at least 48 hours after the first specimen.
- It has been shown that immunological pregnancy tests may yield false results in cases of several medications and diseases such as rheumatoid arthritis. In such cases, the interpretation of the pregnancy test should be done carefully.

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

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It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

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After checking the above mentioned items without finding any error contact your distributor or DRG directly.

## 9 PERFORMANCE CHARACTERISTICS

### 9.1 Assay Dynamic Range

The range of the assay is between 0.44 – 200 mIU/mL.

### 9.2 Specificity of Antibodies (Cross-Reactivity)

The following substances were tested for cross reactivity of the assay:

Protein	Concentration	Produced Colour Intensity Equivalent to β-HCG in serum (mIU/mL)
hLH	300 mIU/mL	6
hLH	200 mIU/mL	< 5
hLH	80 mIU/mL	< 5
TSH	75 µIU/mL	8
TSH	50 µIU/mL	< 5
TSH	25 µIU/mL	< 5
FSH	200 mIU/mL	< 5
FSH	50 mIU/mL	< 5

### 9.3 Sensitivity

The analytical sensitivity was calculated from the mean plus two standard deviations of twenty (20) replicate analyses of the Standard 0 and was found to be 0.44 mIU/mL.

### 9.4 Reproducibility

#### 9.4.1 Intra Assay Variation

The within assay variability is shown below:

Sample	n	Mean (mIU/mL)	CV (%)
1	20	5.9	9.4
2	20	18.6	4.0
3	20	147.9	3.5

#### 9.4.2 Inter Assay Variation

The between assay variability is shown below:

Sample	Mean (mIU/mL)	CV (%)
1	5.6	9.9
2	17.1	7.2
3	140.0	6.9

### 9.5 Recovery

Samples have been spiked by adding β-HCG solutions with known concentrations.

The recovery (%) was calculated by multiplying the ratio of measured and expected values with 100.

	Sample 1	Sample 2	Sample 3
Concentration (mIU/mL)	30.9	69.0	115.2
Average Recovery (%)	94.4	91.3	90.7
Range of Recovery (%)	from	91.3	88.6
	to	99.3	94.2

### 9.6 Linearity

Samples were measured undiluted and in serial dilutions with standard 0. The recovery (%) was calculated by multiplying the ratio of expected and measured values with 100.

	Sample 1	Sample 2	Sample 3
Concentration (mIU/mL)	39.7	75.6	128.4
Average Recovery (%)	102.0	95.5	96.8
Range of Recovery (%)	from	90.9	86.2
	to	115.0	109.3

### 10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

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#### 10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 30 mg/mL) have no influence on the assay results.

#### 10.2 Drug Interferences

Until today no substances (drugs) are known to us, which have an influence to the measurement of β-HCG in a sample.

#### 10.3 High-Dose-Hook Effect

No hook effect was observed in this test up to 15 800 mIU/mL of β-HCG.

### 11 LEGAL ASPECTS

#### 11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

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Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

### 12 REFERENCES / LITERATURE

1. Sancken, U., Bartels, I., Der sogenannte Triple-Test. Reproduktionsmedizin 15:276-284, 1999
2. Cole L.A. Clinical Chemistry 43:12, 2233-2243, 1997
3. Iles R.K., Javid M.K., Gunn L.K., Chard T. Clinical Chemistry 45:4, 532-538, 1999
4. Iles R.K., Javid M.K., Gunn L.K., Chard T. Clinical Chemistry 40:3, 484-485, 1994
5. Korhonen J., Alifhan H., Ylöstalo P., Veldhuis J., Stenman U.-H. Clinical Chemistry 43:11, 2155-2163, 1997

**Free  $\beta$ -Subunit Human Chorionic Gonadotropin (Free Beta hCG) Test System**  
**Product Code: 2025-300**

**1.0 INTRODUCTION**

Intended Use: The Quantitative Determination of Free Beta ( $\beta$ ) Chorionic Gonadotropin Subunit Concentration in Human Serum by a Microplate Enzyme Immunoassay, Colorimetric

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Human chorionic gonadotropin (hCG) concentration increases dramatically in blood and urine during normal pregnancy. hCG is secreted by placental tissue, beginning with the primitive trophoblast, almost from the time of implantation, and serves to support the corpus luteum during the early weeks of pregnancy. hCG or hCG similar glycoproteins can also be produced by a wide variety of trophoblastic and nontrophoblastic tumors. The measurement of hCG, by assay systems with suitable sensitivity and specificity has proven great value in the detection of pregnancy and the potential diagnosis of early pregnancy disorders. Free  $\beta$ -subunit hCG testing has improved the diagnostic probability of abnormal pregnancy/disease states.<sup>1</sup>

Patients with trophoblastic diseases produce ordinary and irregular forms of hCG; e.g. nicked hCG, hCG missing the  $\beta$ -subunit C-terminal segment, hyperglycosylated hCG and free  $\beta$ -subunit. On the other hand, common epithelial tumors of the urogenital tract frequently express the free  $\beta$ -subunit of hCG with no concomitant expression of its heterodimer partner, the common  $\alpha$ -subunit of the glycoprotein hormone. While most hCG assays do a very good job of monitoring the normal pregnancies, still there needs to be a system of differential diagnosis of ovarian tumors, epithelial tumors and trophoblastic malfunctions. That is where determination of free  $\alpha$ -subunit, free  $\beta$ -subunit, nicked hCG and nonnicked hCG etc are of individual value.

Although Free  $\beta$ -Subunit hCG normally constitutes less than 1% of the total hCG concentrations in normal pregnancy, it constitutes a significant part (as much as 26% of hCG) in trophoblastic disease.<sup>2,3</sup> There is also increasing evidence that free beta subunit may be better than total hCG measurement in assessing Down's Syndrome.<sup>4</sup>

In this method, Free Beta hCG calibrator, patient specimen or control is first added to an anti-Free Beta hCG coated well. Assay Buffer is added and the reagents mixed. Reaction between the antibody and native Free Beta hCG forms complex that binds with the antibodies coated to the well. The excess serum proteins are washed away via a wash step. Another enzyme labeled monoclonal antibody specific to Free Beta hCG is added to the wells. The enzyme labeled antibody binds to the Free Beta hCG already immobilized on the well. Excess enzyme is washed off via a wash step. A color is generated by the addition of a substrate. The intensity of the color generation is directly proportional to the concentration of the Free Beta hCG in the sample.


The employment of several serum references of known Free Beta Chorionic Gonadotropin levels permits the construction of a dose response curve of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with Free Beta hCG concentration.

**3.0 PRINCIPLE**

**Immunoenzymometric sequential assay (TYPE 4):**


The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitope recognition, in excess, and native antigen.

Upon mixing assay buffer and a serum containing the native antigen, reaction results between the native antigen and the coated antibody, forming an antibody-antigen complex. This interaction is illustrated below:

$Ab_{(c)}$  = Coated Antibody (Excess Quantity)  
 $Ag$  = Native Antigen (Variable Quantity)  
 $AgAb_{(c)}$  = Antigen-Antibody complex (Variable Quant.)  
 $k_a$  = Rate Constant of Association  
 $k_a$  = Rate Constant of Disassociation

After a suitable incubation period, the antibody-antigen bound fraction is separated from unbound antigen by decantation or aspiration. Another antibody (directed at a different epitope) labeled with an enzyme is added. Another interaction occurs to form an enzyme labeled antibody-antigen-biotinylated-antibody complex on the surface of the wells.

$EnzAb$  = Enzyme labeled Antibody (Excess Quantity)  
 $k_b$  = Rate Constant of Association  
 $k_b$  = Rate Constant of Dissociation

Excess enzyme is washed off via a wash step. A suitable substrate is added to produce color measurable with the use of a microplate spectrophotometer. The enzyme activity on the well is directly proportional to the native free antigen concentration. By utilizing several different serum references of known antigen concentration, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

**4.0 REAGENTS**

**Materials Provided:**

**A. Free Beta hCG Calibrators – 1 ml/vial - Icons A-F**

Six (6) vials of references  $\beta$ -hCG Antigen at levels of 0(A), 10(B), 25(C), 50(D), 100(E) and 250(F) ng/ml. Store at 2-8°C. A preservative has been added.

**Note:** The calibrators, human serum based, were calibrated using a reference preparation, which was assayed against the WHO 1<sup>st</sup> IRP (75/551).

**Conversion to mass units =** One (1) mIU/ml is equivalent to 1 ng/ml

**B. Assay Buffer – 7 ml/vial - Icon B**

One (1) vial containing buffer, dye, surfactants and preservatives. Store at 2-8°C.

**C. Free Beta hCG Enzyme Reagent – 13 ml/vial - Icon C**

One (1) vial containing Enzyme (HRP) labeled monoclonal mouse IgG in buffer, dye, and preservative. Store at 2-8°C.

**D. Free Beta hCG Antibody Coated Plate – 96 wells - Icon D**

One 96-well microplate coated with anti- $\beta$ hCG mouse IgG and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

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**G. Substrate B – 7 ml/vial - Icon S<sup>B</sup>**

One (1) vial containing hydrogen peroxide (HB<sub>2B</sub>OB<sub>2B</sub>) in buffer. Store at 2-8°C.

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**Note 2:** Opened reagents are stable for sixty (60) days when stored at 2-8°C. **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on the label.**

**Note 3:** Above reagents are for a single 96-well microplate. For other kit configurations, refer to table at the end of the instructions.

**4.1 Required But Not Provided:**

- Pipette capable of delivering 0.025ml (25 $\mu$ l), 0.050ml (50 $\mu$ l) 0.100 ml (100  $\mu$ l) volumes with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.50 ml (50  $\mu$ l), 0.100 ml (100  $\mu$ l) and 0.350 ml (0.350  $\mu$ l) volumes with a precision of better than 1.5%.
- Microplate washers or a squeeze bottle (optional).
- Microplate Reader with 450nm and 620nm wavelength absorbance capability.
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplates cover for incubation steps.
- Vacuum aspirator (optional) for wash steps.
- Timer.
- Quality control materials.

**5.0 PRECAUTIONS**

**For In Vitro Diagnostic Use  
 Not for Internal or External Use in Humans or Animals**

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The specimens shall be blood serum in type and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain redtop venipuncture tube without additives or anti-coagulants. Allow the blood to clot for samples. Centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.050 ml (50  $\mu$ l) of the specimen is required.

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**8.0 REAGENT PREPARATION**

- Wash Buffer**  
Dilute contents of wash concentrate to 1000ml with distilled or deionized water in a suitable storage container. Store diluted buffer at 2-30°C for up to 60 days.
- Working Substrate Solution** – Stable for one (1) year  
Pour the contents of the amber vial labeled Solution 'A' into the clear vial labeled Solution 'B'. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2 - 8°C.

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- Format the microplates' wells for each serum reference, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
- Pipette 0.025 ml (25  $\mu$ l) of the appropriate serum reference, control or specimen into the assigned well.
- Add 0.050ml of Assay Buffer to each well. **It is very important to dispense all reagents close to the bottom of the coated well.**
- Swirl the microplate gently for 20-30 seconds to mix and cover.
- Incubate 30 minutes at room temperature.
- Discard the contents of the microplate by decantation or aspiration. If decanting, tap and blot the plate dry with absorbent paper.
- Add 0.350 ml (350  $\mu$ l) of wash buffer (see Reagent Preparation section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100  $\mu$ l) of the Free Beta hCG Enzyme Reagent to each well.  
**DO NOT SHAKE THE PLATE AFTER ENZYME ADDITION**
- Cover and incubate 15 minutes at room temperature.
- Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
- Add 0.350 ml (350  $\mu$ l) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100  $\mu$ l) of working substrate solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells.**  
**DO NOT SHAKE PLATE AFTER SUBSTRATE ADDITION**
- Incubate at room temperature for fifteen (15) minutes.
- Add 0.050 ml (50  $\mu$ l) of stop solution to each well and mix gently for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
- Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. **The results should be read within thirty (30) minutes of adding the stop solution.**  
**Note: For re-assaying specimens with concentrations greater than 250 ng/ml, see 12.1(#9) "Assay Performance" for dilution details.**

**10.0 CALCULATION OF RESULTS**

A dose response curve is used to ascertain the concentration of Free Beta hCG in unknown specimens.

- Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
- Plot the absorbance for each duplicate serum reference versus the corresponding Free Beta hCG concentration in ng/ml on linear graph paper (do not average the duplicates of the serum references before plotting).
- Draw the best-fit curve through the plotted points.
- To determine the concentration of Free Beta hCG for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (0.366) intersects the dose response curve at 14.7 ng/ml Free Beta hCG concentration (See Figure 1).

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
- Maximum Absorbance (Calibrator 'F') = >1.3
- Maximum Absorbance (Calibrator 'A') =  $\leq U \leq 0.1$
- Four out of six quality control pools should be within the established ranges.

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### 12.1 Assay Performance

- It is important that the time of reaction in each well is held constant to achieve reproducible results.
- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- The reagents for the test system procedures have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays (Boscato LM, Stuart MC. "Heterophilic antibodies: a problem for all immunoassays" *Clin.Chem.* 1988:3427-33). For diagnostic purposes, the results from this assay should be used in combination with clinical examination, patient history and all other clinical findings.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Patient specimens with Free Beta hCG concentrations above 250 ng/ml may be diluted (for example 1/10) with normal male serum (Free Beta hCG < 1 ng/ml) and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind's IFU may yield inaccurate results.
- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
- Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from [Monobind@monobind.com](mailto:Monobind@monobind.com).

### 12.2 Interpretation

**1. Measurements and interpretation of results must be performed by a skilled individual or trained professional.**

- Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- The reagents for AccuBind® ELISA procedures have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays (Boscato LM, Stuart MC. "Heterophilic antibodies: a problem for all immunoassays" *Clin. Chem* (1988:3427-33). For diagnostic purposes, the results from this assay should be in combination with clinical examination, patient history, and all other clinical findings.
- For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability.**
- If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

### 13.0 EXPECTED RANGES OF VALUES

Serum Free beta hCG (as intact hCG) increases rapidly in normal pregnancy reaching maximum levels of approximately 60ng/ml

(which is equivalent to 60 mIU/ml) at eight-nine weeks of gestation. This is followed by a gradual decline during the next eleven to twelve weeks. The ration of Free Beta hCG to intact hCG reaches 1% at five weeks of pregnancy and remains constant at approximately 0.5% (w/w) for the remaining twenty-two weeks.<sup>2</sup>

The use of Free Beta hCG in combination with AFP levels as a screening protocol for Down syndrome (Trisomy 21) has been promoted to achieve high detection efficiency with low false positive rates.<sup>5</sup>

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## 14.0 PERFORMANCE CHARACTERISTICS

### 14.1 Precision

The within and between assay precisions of the Free Beta hCG AccuBind® ELISA test system were determined by analyses on three different levels of control sera. The number, mean value, standard deviation and coefficient of variation for each of these control sera are presented in Table 3 and Table 4.

TABLE 2  
Within Assay Precision (Values in ng/ml)

Sample	N	X	σ	C.V.
Level 1	20	3.5	0.15	4.3%
Level 2	20	10.5	0.55	5.2%
Level 3	20	49.6	1.50	3.0%

TABLE 3  
Between Assay Precision\* (Values in ng/ml)

Sample	N	X	σ	C.V.
Level 1	10	3.1	0.17	5.5%
Level 2	10	11.2	0.71	6.3%
Level 3	10	48.5	1.75	3.6%

\*As measured in ten experiments in duplicate.

### 14.2 Sensitivity

The Free Beta hCG AccuBind® ELISA test system has a sensitivity of 0.0007 mIU. This is equivalent to a sample containing 0.026 ng/ml Free Beta hCG concentration.

### 14.3 Specificity

The cross-reactivity of the Free Beta hCG AccuBind® ELISA test method to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of Chorionic Gonadotropin needed to produce the same absorbance.

Substance	Cross Reactivity
fβ-hCG subunit	1.0000
Intact Chorionic Gonadotropin (hCG)	< 0.0001
Follicle Stimulating Hormone (FSH)	< 0.0001
Luteinizing Hormone (LH)	< 0.0001
Thyroid Stimulating Hormone (TSH)	< 0.0001

## 15.0 REFERENCES

- Henry JB, *Clinical Diagnosis and Management by Laboratory Methods*, WB Saunders Company, 486 (1996).
- Osturk M, et al, *Endocrinology*, 120, 549 (1987).
- Cole LA, et al, *Serono Symposia Publ*, 65, 59-78 (1989).
- Macri JN, and Spencer K, *Am J Obstet Gynecol*, 174, 1668-69 (1996).
- Macri JN, et al, *Am J Obstet Gynecol*, 163, 1248 (1990).
- Cole LA, "β core fragment (β-core, UGP or UGF)", *Tumor Marker Update*, 6, 69-75 (1994).
- Yamanaka N, Kawabata G, Morisue K, Hazama M, Nishimura R, "Urinary hCG β-core fragment as a tumor marker for bladder

- cancer", *Nippon Hinyokika Gakkai Zasshi*, 84, 700-706 (1993).
- Kinugasa M, Nishimura R, Hasegawa K, Okamura M, Kimura A, Ohtsu F", Assessment of urinary β-core fragment of hCG as a tumor marker of cervical cancer", *Acta Obstet Gynecol Jpn*, 11, 188-94 (1992).
  - Cole LA, "Stability of hCG free β-subunit in urine", *Prenat Diagnosis*, 17, 185-89 (1996).
  - Sancken U, Bahner D, "The effect of thermal instability of intact human chorionic gonadotropin on the application of its free β-subunit as a serum marker of Down's syndrome screening", *Prenat Diagnosis*, 15, 731-738 (1995).
  - Javadpour N, "Current status of tumor markers in testicular cancer", *Eur Urol*, 21, 34-36 (1992).
  - Canick JA, Kellner DN Jr, Palomaki GE, Walker RP, Osathanondh R, "Second trimester levels of maternal urinary gonadotropin peptide in Down syndrome pregnancy", *Prenat Diagnosis*, 15, 752-759 (1995).

Effective Date: 2024-Jun-21 Rev. 3 DCO: 1678  
MP2025 Product Code: 2025-300

Size	96(A)	192(B)
Reagent (fill)	A) 1ml set	1ml set
	B) 1 (7ml)	2 (7ml)
	C) 1 (13ml)	2 (13ml)
	D) 1 plate	2 plates
	E) 1 (20ml)	1 (20ml)
	F) 1 (7ml)	2 (7ml)
	G) 1 (7ml)	2 (7ml)
	H) 1 (8ml)	2 (8ml)

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
### 1. INTENDED USE

---

The *Aspergillus fumigatus* IgG ELISA is intended for the qualitative determination of IgG class antibodies against *Aspergillus fumigatus* in human serum or plasma (citrate, heparin).

### 2. PRINCIPLE OF THE ASSAY

---

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

### 3. MATERIALS

---

#### 3.1. Reagents supplied

- **Microtiterplate:** 12 break-apart 8-well snap-off strips coated with *Aspergillus fumigatus* antigens; in resealable aluminium foil.
- **DIL:** 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; coloured yellow; ready to use; white cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).
- **SOLN STOP:** 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- **WASH BUF 20x:** 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M), pH 7.2 ± 0.2, for washing the wells; white cap; 0.2% (w/v) 5-Bromo-5-nitro-1,3-dioxane.
- **Conjugate:** 1 bottle containing 20 mL of peroxidase labelled antibody to human IgG in phosphate buffer (10 mM); coloured blue; ready to use; black cap.
- **SUB TMB:** 1 bottle containing 15 mL 3,3',5,5'-tetramethylbenzidine (TMB), < 0.1 %; ready to use; yellow cap.
- **Positive Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; red cap; ≤ 0.02% (v/v) MIT.
- **Cut-off Control:** 1 vial containing 3 mL control; coloured yellow; ready to use; green cap; ≤ 0.02% (v/v) MIT.
- **Negative Control:** 1 vial containing 2 mL control; coloured yellow; ready to use; blue cap; ≤ 0.0015% (v/v) CMIT/ MIT (3:1).

For hazard and precautionary statements see 11.1.

#### 3.2. Materials supplied

- 1 Cover foil
- 1 Instructions for use (IFU)

#### 3.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm
- Incubator 37 °C
- Manual or automatic equipment for rinsing Microtiterplates
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

### 4. STABILITY AND STORAGE

---

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

### 5. REAGENT PREPARATION

---

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

#### 5.1. Microtiterplate

The break-apart snap-off strips are coated with *Aspergillus fumigatus* antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

#### 5.2. WASH BUF 20x

Dilute WASH BUF 20x 1 + 19; e. g. 10 mL WASH BUF 20x + 190 mL distilled water. The diluted buffer (WASH BUF 1x) is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g., in a water bath. Mix well before dilution.

### 5.3. **SUB TMB**

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. **SUB TMB** should be colourless or could have a slight blue tinge. If **SUB TMB** turns into blue, it may have become contaminated and should be thrown away.

## 6. SAMPLE COLLECTION AND PREPARATION

---

Use human serum or plasma (citrate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise, they should be aliquoted and stored deep-frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

### 6.1. Sample Dilution

Before assaying, all samples should be diluted 1+100 with **DIL**. Dispense 10 µL sample and 1 mL **DIL** into tubes to obtain a 1+100 dilution and thoroughly mix with a Vortex.

## 7. ASSAY PROCEDURE

---

Please read the instructions for use carefully **before** performing the assay. Result reliability depends on strict adherence to the instructions for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of **WASH BUF 1x** from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 11. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to 37 ± 1 °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. **Incubate for 60 min at 37 ± 1 °C.**
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of **WASH BUF 1x**. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step!  
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. **Incubate for 30 min at room temperature (20...25 °C).** Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL **SUB TMB** into all wells.
9. **Incubate for exactly 15 min at room temperature (20...25 °C) in the dark.** A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL **SOLN STOP** into all wells in the same order and at the same rate as for **SUB TMB**, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of **SOLN STOP**.

### 7.1. Measurement

Adjust the ELISA Microtiterplate reader to zero using the **Substrate Blank**.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

**Measure the absorbance** of all wells at **450 nm** and record the absorbance values for each standard/control and sample.

Bichromatic measurement using a reference wavelength of 620 nm is recommended.

**Where applicable calculate the mean absorbance values of all duplicates.**

## 8. RESULTS

---

### 8.1. Run Validation Criteria

In order for an assay run to be considered valid, these instructions for use have to be strictly followed and the following criteria must be met:

- **Substrate Blank:** Absorbance value < **0.100**
- **Negative Control:** Absorbance value < **0.200** and < **Cut-off**
- **Cut-off Control:** Absorbance value **0.150 – 1.300**
- **Positive Control:** Absorbance value > **Cut-off**

If these criteria are not met, the test is not valid and must be repeated.

## 8.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 = 0.86 / 2 = 0.43  
Cut-off = 0.43

### 8.2.1. Results in Units [NTU]

$$\frac{\text{Sample (mean) absorbance value} \times 10}{\text{Cut-off}} = [\text{NovaTec Units} = \text{NTU}]$$

Example: 
$$\frac{1.591 \times 10}{0.43} = 37 \text{ NTU}$$

## 8.3. Interpretation of Results

Cut-off	10 NTU	-
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks.
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.  
In immunocompromised patients and newborns serological data only have restricted value.

### 8.3.1. Antibody Isotypes and State of Infection

Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer with low IgG titer: → suggests a current or very recent infection Rare: → persisting IgM
IgG	Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: → may indicate a past infection

## 9. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

### 9.1. Precision

Intraassay	n	Mean (OD)	CV (%)
#1	24	0.408	5.21
#2	24	1.181	3.01
#3	24	0.945	1.52
Interassay	n	Mean (NTU)	CV (%)
#1	12	29.72	6.25
#2	12	22.93	4.67
#3	12	4.93	4.98

### 9.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 95.65% (95% confidence interval: 78.05% - 99.89%).

### 9.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 100% (95% confidence interval: 83.16% - 100%).

### 9.4. Interferences

Interferences with hemolytic, lipemic or icteric samples are not observed up to a concentration of 10 mg/mL hemoglobin, 5 mg/mL triglycerides and 0.5 mg/mL bilirubin.

### 9.5. Cross Reactivity

Investigation of a sample panel with antibody activities to potentially cross-reacting parameters (rheumatoid factors and antibodies to *Candida albicans*) did not reveal evidence of false-positive results due to cross-reactions.

## 10. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.

## 11. PRECAUTIONS AND WARNINGS

---

- The test procedure, the information, the precautions and warnings in the instructions for use have to be strictly followed. The use of the testkits with analyzers and similar equipment has to be validated. Any change in design, composition and test procedure as well as for any use in combination with other products not approved by the manufacturer is not authorized; the user himself is responsible for such changes. The manufacturer is not liable for false results and incidents for these reasons. The manufacturer is not liable for any results by visual analysis of the patient samples.
- Only for in-vitro diagnostic use.
- All materials of human or animal origin should be regarded and handled as potentially infectious.
- All components of human origin used for the production of these reagents have been tested for anti-HIV antibodies, anti-HCV antibodies and HBsAg and have been found to be non-reactive.
- Do not interchange reagents or Microtiterplates of different production lots.
- No reagents of other manufacturers should be used along with reagents of this test kit.
- Do not use reagents after expiry date stated on the label.
- Use only clean pipette tips, dispensers, and lab ware.
- Do not interchange screw caps of reagent vials to avoid cross-contamination.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.
- After first opening and subsequent storage check conjugate and standard/control vials for microbial contamination prior to further use.
- To avoid cross-contamination and falsely elevated results pipette patient samples and dispense reagents without splashing accurately into the wells.
- The ELISA is only designed for qualified personnel following the standards of good laboratory practice (GLP).
- For further internal quality control each laboratory should additionally use known samples.

### 11.1. Safety note for reagents containing hazardous substances

Reagents may contain CMIT/MIT (3:1) or MIT (refer to 3.1).  
Therefore, the following hazard and precautionary statements apply.

	<b>Warning</b>	H317	May cause an allergic skin reaction.
		P261	Avoid breathing spray.
		P280	Wear protective gloves/ protective clothing.
		P302+P352	IF ON SKIN: Wash with plenty of soap and water.
		P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
		P362+P364	Take off contaminated clothing and wash it before reuse.

Reagents may contain 5-Bromo-5-nitro-1,3-dioxane (refer to 3.1).  
Therefore, the following hazard and precautionary statements apply.

	<b>Warning</b>	H315	Causes skin irritation.
		H319	Causes serious eye irritation.
		P280	Wear protective gloves/ protective clothing.
		P302+P352	IF ON SKIN: Wash with plenty of soap and water.
		P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
		P337+P313	If eye irritation persists: Get medical advice/attention.

Further information can be found in the safety data sheet.

### 11.2. Disposal Considerations

Residues of chemicals and preparations are generally considered as hazardous waste. The disposal of this kind of waste is regulated through national and regional laws and regulations. Contact your local authorities or waste management companies which will give advice on how to dispose hazardous waste.

For information about the packaging materials refer to PACKAGING MATERIALS.

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Aspergillus fumigatus IgG

## Test Preparation

Prepare reagents and samples as described.  
Establish the distribution and identification plan for all samples and standards/controls.  
Select the required number of microtiter strips or wells and insert them into the holder.

## Assay Procedure

	Substrate Blank (A1)	Negative Control	Cut-off Control	Positive Control	Sample (diluted 1+100)
Negative Control	-	100 µL	-	-	-
Cut-off Control	-	-	100 µL	-	-
Positive Control	-	-	-	100 µL	-
Sample (diluted 1+100)	-	-	-	-	100 µL
Cover wells with foil supplied in the kit <b>Incubate for 1 h at 37 ± 1 °C</b> Wash each well three times with 300 µL of of <b>WASH   BUF   1x</b>					
Conjugate	-	100 µL	100 µL	100 µL	100 µL
<b>Incubate for 30 min at room temperature (20...25 °C)</b> Do not expose to direct sunlight Wash each well three times with 300 µL of of <b>WASH   BUF   1x</b>					
<b>SUB   TMB</b>	100 µL	100 µL	100 µL	100 µL	100 µL
<b>Incubate for exactly 15 min at room temperature (20...25 °C) in the dark</b>					
<b>SOLN   STOP</b>	100 µL	100 µL	100 µL	100 µL	100 µL
Photometric measurement at 450 nm (reference wavelength: 620 nm)					

### **Gold Standard Diagnostics Frankfurt GmbH**

Waldstrasse 23 A6

63128 Dietzenbach, Germany

Tel.: +49 6074 23698-0

Fax: +49 6074 23698-900

E-Mail: [info.frankfurt@eu.goldstandarddiagnostics.com](mailto:info.frankfurt@eu.goldstandarddiagnostics.com)

Website: [clinical.goldstandarddiagnostics.com](http://clinical.goldstandarddiagnostics.com)

ASPG0680\_IFU\_rev01\_fromLot\_055N

## Unconjugated Estriol (uE3) Test System Product Code: 5025-300

### 1.0 INTRODUCTION

**Intended Use: The Quantitative Determination of Unconjugated (Free) Estriol Concentration in Human Serum or Plasma by a Microplate Enzyme Immunoassay, Colorimetric**

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The last few years have seen the development of screening for fetal Down Syndrome by measurement of multiple markers in maternal blood (1). Although amniocentesis has been widely available for more than 40 years it can only be selectively used to diagnose the disorder because of the hazard to fetus. Of most employed for differential diagnosis the commonly used procedures are AFP, hCG, free beta-HCG and unconjugated estriol (2).

Unconjugated estriol in the serum of pregnant women originates almost exclusively from precursors in the fetus, via the placenta (3). The clinical evidence shows that in uncomplicated pregnancies, the production of estriol increases steadily throughout the last trimester; however, in pregnancies complicated by placental insufficiency the synthesis of estriol decreases rapidly. For many years the most commonly used method for monitoring estriol synthesis (as an index to fetal stress) has been to measure estriol and estriol conjugates in a 24 hr urine sample (4). However, changes in renal clearance and diurnal variations can make the results of these determinations suspect. In recent years investigators have found the determinations of unconjugated estriol in pregnancy plasma, as an alternative to the urinary assay, to be a better marker of fetal stress (5). Abnormally low levels of estriol in a pregnant woman may indicate a problem with the development of the child. Levels of estriol in non-pregnant women do not change much after menopause, and levels are not significantly different from levels in men (6).

The Monobind unconjugated estriol EIA Kit uses a specific anti-estriol antibody, and does not require prior sample extraction of serum or plasma. Cross-reactivity to other naturally-occurring and structurally related steroids is low.


The employment of several serum references of known Estriol concentration permits construction of a graph of activity and concentration. From comparison to the dose response curve, an unknown specimen's activity can be correlated with Estriol concentration.

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
The essential reagents required for an enzyme immunoassay include antibody, enzyme-antigen conjugate and native antigen. Upon mixing biotinylated antibody, enzyme-antigen conjugate and a serum containing the native antigen, a competition reaction results between the native antigen and the enzyme-antigen

conjugate for a limited number of antibody binding sites. The interaction is illustrated by the following equation:

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### 4.0 REAGENTS

#### A. uE3 Calibrators – 1ml/vial - Icons A-F

Six (6) vials of serum reference for unconjugated estriol at concentrations of 0 (A), 0.4 (B), 2.0 (C), 5.0 (D), 15 (E), and 30.0 (F) in ng/ml. Store at 2-8°C. A preservative has been added. The calibrators can be expressed in molar concentrations (nM/L) by the conversion factor 3.45. For example: 1ng/ml x 3.45 = 3.45 nM/L

#### B. uE3 Enzyme Reagent – 6.0 ml/vial – Icon $\text{E}$

One (1) vial contains Estriol (Analog)-horseradish peroxidase (HRP) conjugate in a protein stabilizing matrix with red dye. Store at 2-8°C.

#### C. uE3 Biotin Reagent – 6.0 ml - Icon $\nabla$

One (1) vial contains anti-unconjugated Estriol biotinylated purified rabbit IgG conjugate in buffer, blue dye and preservative. Store at 2-8°C.

#### D. Streptavidin Coated Plate – 96 wells –Icon $\downarrow$

One 96-well microplate coated with 1.0  $\mu\text{g/ml}$  streptavidin and packaged in an aluminum bag with a drying agent. Store at 2-8°C.

#### E. Wash Solution Concentrate – 20ml/vial - Icon $\text{D}$

One (1) vial contains a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.

#### F. Substrate A – 7ml/vial - Icon $\text{S}^{\text{A}}$

One (1) vial contains tetramethylbenzidine (TMB) in buffer. Store at 2-8°C.

#### G. Substrate B – 7ml/vial - Icon $\text{S}^{\text{B}}$

One (1) vial contains hydrogen peroxide ( $\text{H}_2\text{O}_2$ ) in buffer. Store at 2-8°C.

#### H. Stop Solution – 8ml/vial - Icon $\text{STOP}$

One (1) vial contains a strong acid (1N HCl). Store at 2-8°C.

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#### 4.1 Required But Not Provided:

- Pipette capable of delivering 0.025ml (25 $\mu\text{l}$ ) and 0.050ml (50 $\mu\text{l}$ ) with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.100ml (100 $\mu\text{l}$ ) and 0.350ml (350 $\mu\text{l}$ ) volumes with a precision of better than 1.5%.
- Adjustable volume (200-1000 $\mu\text{l}$ ) dispenser(s) for conjugate.
- Microplate washer or a squeeze bottle (optional).
- Microplate Reader with 450nm and 620nm wavelength absorbance capability.
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplate cover for incubation steps.

8. Vacuum aspirator (optional) for wash steps.

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The specimens shall be blood; serum or heparinized plasma in type and taken with the usual precautions in the collection of venipuncture samples. The blood should be collected in a redtop (with or without gel additives) venipuncture tube or for plasma use evacuated tube(s) containing heparin. Allow the blood to clot for serum samples. Centrifuge the specimen to separate the serum or plasma from the cells.

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Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.050ml (50 $\mu\text{l}$ ) of the specimen is required.

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### 8.0 REAGENT PREPARATION

- Wash Buffer**  
Dilute contents of wash solution to 1000ml with distilled or deionized water in a suitable storage container. Diluted buffer can be stored at 2-30°C for up to 60 days.
- Working Substrate Solution - Stable for 1 year**  
Pour the contents of the amber vial labeled Solution 'A' into the clear vial labeled Solution 'B'. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2 - 8°C.

**Note1 : Do not use the working substrate if it looks blue.**

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- Format the microplates' wells for each calibrator, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
- Pipette 0.025 ml (25 $\mu\text{l}$ ) of the appropriate calibrator, control or specimen into the assigned well.
- Add 0.050 ml (50 $\mu\text{l}$ ) of the uE3 Enzyme Reagent to all wells (see Reagent Preparation Section).
- Swirl the microplate gently for 20-30 seconds to mix.
- Add 0.050 ml (50 $\mu\text{l}$ ) of uE3 Biotin Reagent to all wells.
- Swirl the microplate gently for 20-30 seconds to mix.
- Cover and incubate for 60 minutes at room temperature.
- Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
- Add 0.350ml (350 $\mu\text{l}$ ) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- Add 0.100 ml (100 $\mu\text{l}$ ) of working substrate solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells.**
- DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION**
- Incubate at room temperature for fifteen (15) minutes.
- Add 0.050ml (50 $\mu\text{l}$ ) of stop solution to each well and gently mix for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
- Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm). **The results should be read within thirty (30) minutes of adding the stop solution.**

**Note:** Dilute the sample, suspected of concentrations higher than 30ng/ml, by diluting 1:2 and/or 1:5 with unconjugated estriol '0' ng/ml calibrator or male patient sera with a known low value for estriol. Multiply the result by the dilution factor of 2 or 5 as required to obtain the concentration of the sample.

### 10.0 CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of unconjugated estriol in unknown specimens.


- Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
- Plot the absorbance for each duplicate serum reference versus the corresponding unconjugated estriol concentration in ng/ml on linear graph paper.
- Connect the points with a best-fit curve.
- To determine the concentration of unconjugated estriol for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance (0.967) intersects the dose response curve at (4.71 ng/ml) unconjugated estriol concentration (See Figure 1).

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**12.1 Assay Performance**

- It is important that the time of reaction in each well is held constant to achieve reproducible results.
- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Patient specimens with Unconjugated E3 concentrations above 30 ng/ml may be diluted (1/2, 1/5 or higher) with Unconjugated E3 "0" calibrator and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind IFU may yield inaccurate results.
- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must

be strictly followed to ensure compliance and proper device usage.

- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.

- Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from [Monobind@monobind.com](mailto:Monobind@monobind.com).

**12.2 Interpretation**

- Measurements and interpretation of results must be performed by a skilled individual or trained professional.**
- Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- The reagents for the test system have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays (Boscato LM, Stuart MC. 'Heterophilic antibodies: a problem for all immunoassays' Clin. Chem. 1988;34:27-33). For diagnostic purposes, the results from this assay should be in combination with clinical examination, patient history and all other clinical findings.
- For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability.**
- If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

**13.0 EXPECTED RANGES OF VALUES**

In agreement with established reference intervals<sup>5</sup> for a "normal" adult population, the expected ranges for the Unconjugated Estriol AccuBind® ELISA Test System are detailed in Table 1.

**TABLE 1**  
**Expected Values for the uE3 EIA Test System (ng/ml)**  
**Male & Non-Pregnant Female < 1.0 ng/ml**

During pregnancy the Unconjugated E3 serum levels rise rapidly till the end of third trimester. (See Table 2 from published Literature).<sup>6</sup>

TABLE 2				
Gestation Week	Expected Range (ng/ml)	Gestation Week	Expected Range (ng/ml)	Twin Pregnancy (ng/ml)
12	0.3-1.0	22	2.7-16.0	3.0-18.0
14	0.4-1.6	26	3.0-18.0	4.0-21.0
16	1.4-6.5	32	4.6-23.0	5.0-25.0
18	1.6-8.5	36	7.2-29.0	7.0-39.0
20	2.1-13.0	40	8.0-39.0	13.0-40.0

It is important to keep in mind that establishment of a range of values which can be expected to be found by a given method for a population of "normal" persons is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the Manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

**14.0 PERFORMANCE CHARACTERISTICS**

**14.1 Precision**

The within and between assay precision of the uE3 AccuBind® ELISA Test System were determined by analyses on three different levels of pool control sera. The number, mean values, standard deviation and coefficient of variation for each of these control sera are presented in Table 3 and Table 4.

**TABLE 3**

Within Assay Precision (Values in ng/ml)				
Sample	N	X	$\sigma$	C.V.
Low	24	1.58	0.13	8.3%
Normal	24	5.17	0.37	7.1%
High	24	9.06	0.59	6.5%

**TABLE 4**

Between Assay Precision (Values in ng/ml)				
Sample	N	X	$\sigma$	C.V.
Low	10	1.47	0.14	9.5%
Normal	10	4.93	0.39	7.9%
High	10	8.99	0.54	6.0%

\*As measured in ten experiments in duplicate over a ten day period.

**14.2 Sensitivity**

The uE3 AccuBind® ELISA Test System has a sensitivity of 2.9 pg/T. This is equivalent to a sample containing a concentration of 0.115 ng/ml. The sensitivity was ascertained by determining the variability of the 0 ng/ml serum calibrator and using the  $2\sigma$  (95% certainty) statistic to calculate the minimum dose.

**14.3 Accuracy**

The uE3 AccuBind® ELISA Test System was compared with a reference method. Biological specimens from low, normal and high Unconjugated Estriol level populations were used; the values ranged from 0.15 – 29.1 ng/ml. The total number of specimens was 158. The least square regression equation and the correlation coefficient were computed for this uE3 in comparison with the reference method. The data obtained is displayed in Table 5.

TABLE 5			
Method	Mean (x)	Least Square Regression Analysis	Correlation Coefficient
Monobind (y)	3.84	y= -0.174+0.979(x)	0.952
Reference (X)	3.74		

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**14.4 Specificity**

The % cross reactivity of the Estriol antibody to selected substances, for determination of Unconjugated Estriol, was evaluated by adding the interfering substance to a serum matrix at massive concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of uE3 needed to displace the same amount of labeled analog.

Substance	% Cross Reactivity
Estriol	100.0000
Androstenedione	0.0001
Cortisol	<0.0001
Cortisone	<0.0001
Corticosterone	<0.0001
DHEA-S	<0.0001
Dihydrotestosterone	0.0001
Estradiol	0.0040
Estriol Glucuronide	<0.0001
Estriol Sulfate	0.6200
Estrone	0.0004
Prednisone	<0.0001
Progesterone	<0.0001
Spirolactone	<0.0001
Testosterone	<0.0001

**15.0 REFERENCES**

- Goebelsman, U, Katagiri, H, Stanczyk et al., Estriol assays in obstetrics. *J. Steroid Biochemistry*, 6, 703-709 (1975).
- Reynolds T, Penny M, "The mathematical basis of multivariate risk screening, with special reference to screening for Down Syndrome associated pregnancy", *Ann Clin Biochem*, 27:452-458 (1990).
- Kastagiri, H., Stanczyk, F and Goebelsman, U, "Estriol in pregnancy.III Development, comparison and use of a specific antisera for rapid radioimmunoassay of unconjugated estriol in pregnancy plasma", *Steroids*, 24, 225 (1974).

- Brown CH, Saffan, BD, Howard C.M and Preedy JR. "The renal clearance of endogenous estrogens in late pregnancy", *J Clinical Investigation*, 43,295 (1964).
- Tietz, NW, ED: *Clinical Guide to Laboratory Tests*, 3<sup>rd</sup> ed. Philadelphia, WA Saunders Co, 1995.
- Cohen M, and Cohen H, "A radioimmunoassay for plasma unconjugated estrogens in normal pregnancy", *Am J. Obstet. Gynecology*, 118, 200 (1974).
- Wright JV, Schliesman B, Robinson L, "Comparative measurements of serum estriol, estradiol and estrone in non-pregnant, postmenopausal women; a preliminary investigation.", *Altern Med Rev*, 4 266-70 (1999).
- NIH State-of-the-Science Conference Statement on Management of Menopause-Related Symptoms. NIH Consensus State Sci Statements. 2005. Mar 21-23; 22(1)1-38.

Rev 4 Date: 2019-Jul-16 DCO: 1303  
MP5025 Product Code: 5025-355

Size	96(A)		192(B)	
	A)	1ml set	B)	1ml set
Reagent (fill)	B)	1 (6ml)	2 (6ml)	
	C)	1 (6ml)	2 (6ml)	
	D)	1 plate	2 plates	
	E)	1 (20ml)	1 (20ml)	
	F)	1 (7ml)	2 (7ml)	
	G)	1 (7ml)	2 (7ml)	
	H)	1 (8ml)	2 (8ml)	

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**Glossary of Symbols**  
(EN 980/ISO 15223)

In Vitro - Diagnostic Medical Device  
 Temperature Limitation Storage Condition (2-8°C)  
 Consult Instructions for Use  
 Catalogue Number  
 Contains Sufficient Test for Σ  
 Batch Code  
 Used By (Expiration Day)  
 Date of Manufacturer  
 Manufacturer

Authorized Rep in European Country  
 European Conformity

**Pregnancy Associated Plasma Protein-A (PAPP-A) Test System**  
**Product Code: 12625-300**

**1.0 INTRODUCTION**

**Intended Use: The Quantitative Determination of PAPP-A Concentration in Human Serum by a Microplate Enzyme Immunoassay, Colorimetric**

**2.0 SUMMARY AND EXPLANATION OF THE TEST**

Trisomy 21, 18, and 13 are three conditions that can cause ill effects to a fetus and, therefore, an infant. Screening for markers and inherent trends has been of interest over the years. One of the proteins of interest for these trisomy disorders is PAPP-A, pregnancy associated plasma protein A. This protein can be found in the serum of every individual, but the concentration rises in maternal serum as gestation time progresses. In studies, a correlation has been shown between the decreased level of PAPP-A and the occurrence of trisomy disorders, particularly Trisomy 21, also known as Down's Syndrome. Along with a few other markers like AFP, uE3, and hCG, trends in PAPP-A have been found to elude to trisomy disorders.<sup>1,2,3,4,5,6</sup>

PAPP-A is produced primarily by the placenta during pregnancy. This glycoprotein has a molecular weight of 740,000 and tends to exist as a heterotetrameric dimer with PromBP, proform major basic protein. The concentration of PAPP-A in maternal blood increases over the time of the pregnancy as the placenta and fetus grows because it is a product of the trophoblast. In general, the concentration level of this protein in maternal serum is indicative of threatened abortion, preterm birth, intrauterine growth restriction, ectopic gravidity, preeclampsia or diabetes mellitus. When tested during the first trimester of pregnancy, PAPP-A is the major marker for Down's syndrome.<sup>3,7,8,9</sup>

The Monobind PAPP-A test system is designed specifically for the testing of the heterotetrameric form important during pregnancy. Another form of PAPP-A also exists in serum, but it is a dimeric form associated with coronary and cardiac conditions. Tests developed for use on pregnancy patients are not designed to test for this dimeric form. When evaluation of this protein is done, it is often compared to the Multiple of Medians (MoM) and represented as a percent of an in-house establish median. Without easy access to an IRP for PAPP-A, this method allows an easier way to compare results between laboratories and testing methods, allowing for a type of reference value to be established.<sup>2,5,8,9</sup>


**3.0 PRINCIPLE**

**Immunoenzymometric Sequential Assay (TYPE 4):**

The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme and immobilized), with different and distinct epitope recognition, in excess, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on the well and


exogenously added biotinylated monoclonal anti- PAPP-A antibody.

Upon mixing monoclonal biotinylated antibody, and a serum containing the native antigen, reaction results between the native antigen and the antibody, forming an antibody-antigen complex. Simultaneously the biotin attached to the antibody binds to the streptavidin coated on the microwells resulting in immobilization of the complex. The interaction is illustrated by the following equation:

$Ag_{(PAPP-A)} - {}^{Biotin}Ab_{(m)} + \text{Streptavidin}_{C.W.} \rightarrow \text{immobilized complex (IC)}$   
 $\text{Streptavidin}_{C.W.}$  = Streptavidin immobilized on well  
 $\text{Immobilized complex (IC)}$  = Ag-Ab bound to the well

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**Materials Provided:**

**A. PAPP-A Calibrators – 0.5 ml/vial – Icons A-F**

Six (6) vials of serum reference for PAPP-A at concentrations of 0 (A), 0.64 (B), 1.6 (C), 3.2 (D), 9.6 (E) and 32 (F) in µg/ml. A preservative has been added. Store at 2-8°C.  
**Note:** To convert µg/ml to mIU/L, multiply by 156.25  
**Ex.** 32 µg/ml x 156.25 = 5000 mIU/L

**B. PAPP-A Enzyme Reagent– 12 ml/vial-Icon<sup>ⓔ</sup>**

One (1) vial of PAPP-A Antibody-horseradish peroxidases (HRP) conjugate in a protein-stabilizing matrix with a preservative. Store at 2-8°C.

**C. PAPP-A Biotin Reagent – 12 ml/vial – Icon ∇**

One (1) vial containing biotinylated monoclonal mouse IgG in buffer, dye, and preservative. Store at 2-8°C.

**Note: Patients must be fasting from biotin therapy/supplements for 8 hours to prevent interference in the test**

**D. PAPP-A Control – 0.5 ml/vial – Icon M**

One (1) vials of serum reference for PAPP-A at concentration established (exact value listed on label). A preservative has been added. Store at 2-8°C.

**E. PAPP-A Diluent – 5.0 ml/vial – Icon U**

One (1) vial of human serum based buffer with salts, surfactants, and preservatives. Store at 2-8°C.

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**G. Wash Solution Concentrate – 20ml/vial – Icon ⬇**

One (1) vial contains a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.

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**I. Stop Solution – 8ml/vial – Icon ⊕**

One (1) vial contains of a strong acid (0.5M H<sub>2</sub>SO<sub>4</sub>). Store at 2-8°C.

**J. Product Instructions.**

**Note 1:** Do not use reagents beyond the kit expiration date.  
**Note 2:** Avoid extended exposure to heat and light. **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on label.**  
**Note 3:** Above reagents are for a single 96-well microplate. For other kit configurations, refer to table at the end of the instructions.

**4.1 Required But Not Provided:**

- Pipette capable of delivering 0.010ml (10µl), 0.050ml (50µl) and 0.100ml (100µl) with a precision of better than 1.5%.
- Dispenser(s) for repetitive deliveries of 0.100ml (100µl) and 0.350ml (350µl) volumes with a precision of better than 1.5%.
- Microplate washer or a squeeze bottle (optional).
- Microplate Reader with 450nm and 620nm wavelength absorbance capability.
- Absorbent Paper for blotting the microplate wells.
- Plastic wrap or microplate covers for incubation steps.
- Vacuum aspirator (optional) for wash steps.
- Timer.
- Quality control materials.

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**6.0 SPECIMEN COLLECTION AND PREPARATION**

The specimens shall be blood, serum in type, and the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained. The blood should be collected in a plain redtop venipuncture tube without additives or anti-coagulants; a red top gel separator tube (also known as a "Tiger" top tube) can be used. Allow the blood to clot. Centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.020ml of the specimen is required.

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**7.0 QUALITY CONTROL**

Each laboratory should assay controls at levels in the low, normal and high range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay

performance limits. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

**8.0 REAGENT PREPARATION**

**1. Wash Buffer**

Dilute contents of wash solution to 1000ml with distilled or deionized water in a suitable storage container. Diluted buffer can be stored at room temperature (2-30°C) for up to 60 days.

**Note 1: Do not use the substrate if it looks blue.**  
**Note 2: Do not use reagents that are contaminated or have bacteria growth.**

**9.0 TEST PROCEDURE**

*Before proceeding with the assay, bring all reagents, serum references and controls to room temperature (20 - 27°C).*

**\*\*Test procedure should be performed by a skilled individual or trained professional\*\***

- Format the microplates' wells for each serum reference, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
- Pipette 0.010ml (10µL) of the appropriate serum reference, control or specimen into the assigned well.
- Add 0.100ml (100µl) of the PAPP-A Biotin Reagent to all wells.
- Swirl the microplate gently for 20-30 seconds to mix.
- Cover and incubate for 30 minutes at room temperature.
- Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
- Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat four (4) additional times for a total of five (5) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- Add 0.100ml (100µl) of PAPP-A Enzyme Reagent to all wells.
- Cover and incubate for 30 minutes at room temperature
- Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
- Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate.
- Repeat four (4) additional times for a total of five (5) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat four (4) additional times.**
- Add 0.100ml (100µl) of substrate solution to all wells. **Always add reagents in the same order to minimize reaction time differences between wells. DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION**
- Incubate at room temperature for 15-20 minutes.
- Add 0.050ml (50µl) of stop solution to each well and gently mix for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
- Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm). **The results should be read within fifteen (15) minutes of adding the stop solution.**

**Note:** Dilute the samples suspected of concentrations higher than 32 µg/ml with PAPP-A diluent and multiply result by dilution factor. For a 1:5 dilution add 40µl of diluent to 10µl of high concentration sample.

## 10.0 CALCULATION OF RESULTS

A dose response curve is used to ascertain the concentration of PAPP-A in unknown specimens.


- Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
- Plot the absorbance for each duplicate serum reference versus the corresponding PAPP-A concentration in µg/ml on linear graph paper (do not average the duplicates of the serum references before plotting).
- Connect the points with a best-fit curve.
- To determine the concentration of PAPP-A for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in µg/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance in the patient/ sample (1.251) intersects the dose response curve at (7.7 µg/ml) PAPP-A concentration (See Figure 1).

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FIGURE 1

\*The above data and table below is for example only. Do not use it for calculating your results.

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- The absorbance (OD) of Calibrator F (32 µg/ml) should be ≥ 1.3.
- Four out of six quality control pools should be within the established ranges.

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- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. Therefore, the substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from Monobind's IFU may yield inaccurate results.
- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be strictly followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
- Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from [Monobind@monobind.com](mailto:Monobind@monobind.com).

### 12.2 Interpretation

- Measurements and interpretation of results must be performed by a skilled individual or trained professional.** Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- The reagents for PAPP-A AccuBind® ELISA have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of Immunoassays (Boscato LM, Stuart MC "Heterophilic antibodies: a problem for all immunoassays" *Clin. Chem* (1988:3427-33). For diagnostic purposes, the results from this assay should be in combination with clinical examination, patient history, and all other clinical findings.
- For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability.**
- If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

## 13.0 EXPECTED RANGES OF VALUES

In agreement with established reference intervals for a "normal" adult population, the expected ranges for the PAPP-A AccuBind® ELISA Test System are detailed in Table 1.

It is recommended to compare values based on the Multiple of Median (MoM) established for the laboratory when assessing patient samples. By dividing the value of the patient sample by the MoM will give a percent value that is used frequently for evaluation.

Expected Values for the PAPP-A Test System	
Term of Gestation (full weeks)	PAPP-A Concentration (µg/ml)
9	3.86
10	7.1
11	10.1
12	16.68
13	23.2

It is important to keep in mind that establishment of a range of values, which can be expected to be found by a given method for a population of "normal" persons, is dependent upon a multiplicity of factors: the specificity of the method, the population tested and the precision of the method in the hands of the analyst. For these reasons each laboratory should depend upon the range of expected values established by the manufacturer only until an in-house range can be determined by the analysts using the method with a population indigenous to the area in which the laboratory is located.

## 14.0 PERFORMANCE CHARACTERISTICS

### 14.1 Precision

The total precision of the PAPP-A AccuBind® Microplate EIA Test System was determined by analyses on six different levels of pool control and patient sera. The mean values, standard deviation and coefficient of variation for each of these control sera are presented in Table 2.

TABLE 2  
Precision data for the PAPP-A Test System

Sample	Mean Value (µg/ml)	Within-Run Precision		Total Precision (n=80)	
		SD	CV%	SD	CV%
		Control 1	0.69	0.03	3.65
Control 2	2.17	0.14	6.29	0.20	9.33
Control 3	8.25	0.42	5.11	0.60	7.29
Patient 1	1.35	0.05	3.68	0.13	9.51
Patient 2	5.12	0.22	4.24	0.43	8.49
Patient 3	28.02	1.51	5.39	2.60	9.29

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### 14.2 Sensitivity

The LoB/LoD/LoQ of the PAPP-A AccuBind® ELISA Test System were calculated in conformance with CLSI EP17-A2 Protocols for Determination of Limits of Detection of Limits of Quantitation. The LoB is 0.019 µg/ml and the LoD=LoQ=0.074 µg/ml.

### 14.3 Accuracy

#### 14.3.1 Linearity

The linearity of the PAPP-A AccuBind® ELISA test system was tested by diluting human serum samples containing high levels of PAPP-A (10 to 37 µg/ml) with low PAPP-A diluent. The system produces excellent linearity through the range of the test up to 37 µg/ml.

#### 14.3.2 Recovery

The recovery of the PAPP-A AccuBind® Microplate ELISA Test System was calculated for five patient samples spiked with 0.5, 2.0, 4.0, 10.0, and 30.0 µg/ml PAPP-A. Recoveries were determined to be within 15% of the expected values for all samples.

#### 14.3.3 Method Comparison

The PAPP-A AccuBind® ELISA Test System was compared with a reference immunoassay method. Biological specimens from low, normal and relatively high PAPP-A level populations were used; the values ranged from 0.1 µg/ml – 36 µg/ml. The total number of such specimens was 50. The least square regression equation and the correlation coefficient were computed for this PAPP-A ELISA in comparison with the reference method.

TABLE 4

Method	Mean	Least Square Regression Analysis	Correlation Coefficient
Monobind (y)	6.52	y = 0.34 + 0.97(x)	0.992
Reference (x)	6.68		

### 14.4 Specificity

The % cross reactivity of the PAPP-A antibody to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between doses of interfering substance to dose of PAPP-A needed to displace the same amount of labeled analog.

TABLE 5

Substance	Cross Reactivity
Free βhCG	ND
AFP	ND
hPRL	ND
FSH	ND

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## Glossary of Symbols (EN 980/ISO 15223)

In Vitro - Diagnostic Medical Device

Temperature Limitation Storage Condition (2-8° C)

Consult Instructions for Use

Catalogue Number

Contains Sufficient Test for Σ

Batch Code

Used By (Expiration Day)

Date of Manufacturer

Manufacturer

**Free Testosterone AccuBind® ELISA Test System**  
**Product Code: 5325-300**  
**Rx Only**

**1.0 INTRODUCTION**

Indications for Use: The device is an Enzyme Immunoassay (EIA) for the quantitative measurement of free testosterone in human serum. Measurement of free testosterone is used in the diagnosis and treatment of disorders involving the male sex hormones (androgens), including primary and secondary hypogonadism, impotence in males and in females; hirsutism (excessive hair) and virilization (masculinization) due to tumors, polycystic ovaries and androgenital syndromes.

**2.0 SUMMARY AND EXPLANATION OF THE TEST**

Testosterone, (17β-Hydroxy-4-androstene-3-one), a C<sub>19</sub> steroid, is the most potent naturally secreted androgen. In normal post pubertal males, testosterone is secreted primarily by the testes, with only a small amount derived from peripheral conversion of 4-Androstene-3, 17-dione (ASD).<sup>2</sup> In adult women, it has been estimated that over 50% of serum testosterone is derived from peripheral conversion of ASD secreted by the adrenal and ovary, with the remainder from direct secretion of testosterone by these glands.

In the male, testosterone is mainly synthesized in the interstitial Leydig cells and the testis, and is regulated by the interstitial cell stimulating hormone (ICSH), or luteinizing hormone (LH) of the anterior pituitary (the female equivalent of ICSH).<sup>3</sup> Testosterone is responsible for the development of secondary sex characteristics, such as the accessory sex organs, the prostate, seminal vesicles and the growth of facial, pubic and axillary hair. Testosterone measurements have been very helpful in evaluating hypogonadal states. Increased testosterone levels in males can be found in complete androgen resistance (testicular feminization). Common causes of decreased testosterone levels in males include: hypogonadism, orchidectomy, estrogen therapy, Klinefelter's syndrome, hypopituitarism, and hepatic cirrhosis.<sup>2,4</sup>

In the female, testosterone levels are normally found to be much lower than those encountered in the healthy male. Testosterone in the female comes from three sources. It is secreted in small quantities by both the adrenal glands and the ovaries, and in healthy women, 50-60% of the daily testosterone production arises from peripheral metabolism of prohormone, chiefly androstenedione. Common causes of increased serum testosterone levels in females include polycystic ovaries (Stein-Leventhal syndrome), ovarian tumors, adrenal tumors and adrenal hyperplasia. Virilization in women is associated with the administration of androgens and endogenous overproduction of testosterone. There appears to be a correlation between serum testosterone levels and the degree of virilization in women, although approximately 25% of women with varying degrees of virilism have serum testosterone levels that fall within the female reference range.


The majority of testosterone is bound to transport proteins: weakly bound to albumin and cortisol binding protein (25-65% females; 45-85% males) and tightly bound to sex hormone-binding globulin (SHBG) (35-75% females; 14-50% males).<sup>5</sup> A small fraction exist as unbound or free testosterone; however, this form is biologically active. Therefore, the free hormone concentration is a better indicator of biological activity than total testosterone.

**3.0 PRINCIPLE**

**Competitive Enzyme Immunoassay (TYPE 5)**

The essential reagents required for a solid phase enzyme immunoassay include immobilized antibody, enzyme-antigen conjugate and native antigen.

Upon mixing immobilized antibody, enzyme-antigen conjugate and a serum containing the native antigen, a competition reaction results between the native antigen and the enzyme-antigen conjugate for a limited number of insolubilized binding sites. The interaction is illustrated by the equation in the following below.

Ab<sub>C.W.</sub> = Monospecific Immobilized Antibody (Constant Quantity)  
 Ag = Native Antigen (Variable Quantity)

EnzAg = Enzyme-antigen Conjugate (Constant Quantity)  
 AgAb<sub>C.W.</sub> = Antigen-Antibody Complex  
 EnzAgAb<sub>C.W.</sub> = Enzyme-antigen Conjugate -Antibody Complex  
 k<sub>a</sub> = Rate Constant of Association  
 k<sub>a</sub> = Rate Constant of Disassociation  
 K = k<sub>a</sub> / k<sub>a</sub> = Equilibrium Constant

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is inversely proportional to the native antigen concentration. By utilizing several different serum references of known antigen concentration, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

**4.0 REAGENTS**

**Materials Provided:**

- A. Free Testosterone Calibrators – 1ml/vial – Icons A-G**  
 Seven (7) vials of serum reference for Free Testosterone at concentrations of 0 (A), 0.2 (B), 1.0 (C), 2.5 (D), 7.5 (E), 20 (F) and 60 (G) in pg/ml. Store at 2-8°C. A preservative has been added. The calibrators can be expressed in molar concentrations (pM/L) by multiplying by 3.47.  
 For example: 1pg/ml x 3.47 = 3.47 pM/L
- B. Free Testosterone Controls – 1ml/vial – Icons L, M, N**  
 Three (3) vials of serum reference for Free Testosterone at low, middle, and high established concentrations (range values listed on labels). A preservative has been added. Store at 2-8 °C.
- C. Free Testosterone Enzyme Reagent – 13ml/vial – Icon E**  
 One (1) vial of Testosterone (Analog)-horseradish peroxidases (HRP) conjugate in a protein stabilizing matrix with dye. Store at 2-8°C.
- D. Free Testosterone Coated Plate – 96 wells – Icon Y**  
 One 96-well microplate coated with testosterone antibody and packaged in an aluminum bag with a drying agent. Store at 2-8°C.
- E. Wash Solution Concentrate – 20ml/vial – Icon D**  
 One (1) vial containing a surfactant in buffered saline. A preservative has been added. Store at 2-8°C.
- F. Substrate A – 7ml/vial – Icon S<sup>A</sup>**  
 One (1) vial contains tetramethylbenzidine (TMB) in buffer. Store at 2-8°C. See "Reagent Preparation."
- G. Substrate B – 7ml/vial – Icon S<sup>B</sup>**  
 One (1) vial contains hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) in buffer. Store at 2-8°C. See "Reagent Preparation."
- H. Stop Solution - 8ml/vial – Icon H**  
 One (1) vial contains a strong acid (1N HCl). Store at 2-8°C.

**I. Product Instructions**

- Note 1:** Do not use reagents beyond the kit expiration date.
- Note 2:** Avoid extended exposure to heat and light. **Opened reagents are stable for sixty (60) days when stored at 2-8°C. Kit and component stability are identified on label.**
- Note 3:** Above reagents are for a single 96-well microplate.

**4.1 Required But Not Provided:**

- 1. Pipette capable of delivering 0.020 & 0.050ml (20µl & 50µl) volumes with a precision of better than 1.5%.
- 2. Dispenser(s) for repetitive deliveries of 1.00 & 0.350ml (100 & 350µl) volumes with a precision of better than 1.5%.
- 3. Adjustable volume (200-1000µl) dispenser(s) for conjugate.
- 4. Microplate washer or a squeeze bottle (optional).
- 5. Microplate Reader with 450nm and 620nm wavelength absorbance capability.
- 6. Absorbent Paper for blotting the microplate wells.
- 7. Plastic wrap or microplate cover for incubation steps.
- 8. Vacuum aspirator (optional) for wash steps.
- 9. Timer.
- 10. Quality control materials.

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All products that contain human serum have been found to be non-reactive for Hepatitis B Surface Antigen, HIV 1&2 and HCV Antibodies by FDA approved tests. Since no known test can offer complete assurance that infectious agents are absent, all human serum products should be handled as potentially hazardous and capable of transmitting disease. Good laboratory procedures for handling blood products can be found in the Center for Disease Control / National Institute of Health, "Biosafety in Microbiological and Biomedical Laboratories," 2nd Edition, 1988, HHS Publication No. (CDC) 88-8395.

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**6.0 SPECIMEN COLLECTION AND PREPARATION**

The specimens shall be blood serum in type, and the usual precautions in the collection of venipuncture samples should be observed. The blood should be collected in a plain redtop venipuncture tube. Allow the blood to clot for serum samples. Centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of five (5) days. If the specimen(s) cannot be assayed within this time, the sample(s) may be stored at temperatures of -20°C for up to 30 days. Avoid use of contaminated devices. Avoid repetitive freezing and thawing. When assayed in duplicate, 0.040ml (40µl) of the specimen is required.

**7.0 QUALITY CONTROL**

Each laboratory should assay controls at levels in the low, normal and high range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits in accordance with local, state, and federal quality control regulations. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

**8.0 REAGENT PREPARATION**

- 1. Wash Buffer**  
 Dilute contents of wash solution to 1000ml with distilled or deionized water in a suitable storage container. Diluted buffer can be stored at 2-30°C for up to 60 days.
- 2. Working Substrate Solution - Stable for 1 year.**  
 Pour the contents of the amber vial labeled Solution 'A' into the clear vial labeled Solution 'B'. Place the yellow cap on the clear vial for easy identification. Mix and label accordingly. Store at 2 - 8°C.

**Note 1: Do not use reagents that are contaminated or have bacteria growth.**  
**Note 2: Do not use the working substrate if it looks blue.**

**9.0 TEST PROCEDURE**

*Before proceeding with the assay, bring all reagents, serum reference calibrators and controls to room temperature (20- 27 °C). \*\*Test Procedure should be performed by a skilled individual or trained professional\*\**


- 1. Format the microplates' wells for each serum reference, control and patient specimen to be assayed in duplicate. **Replace any unused microwell strips back into the aluminum bag, seal and store at 2-8°C.**
- 2. Pipette 0.020ml (20µL) of the appropriate serum reference, control or specimen into the assigned well.
- 3. Add 0.100ml (100µl) of the Free Testosterone Enzyme Reagent to all wells.
- 4. Swirl the microplate gently for 20-30 seconds to mix.
- 5. Cover and incubate for 60 minutes at room temperature.
- 6. Discard the contents of the microplate by decantation or aspiration. If decanting, blot the plate dry with absorbent paper.
- 7. Add 0.350ml (350µl) of wash buffer (see Reagent Preparation Section), decant (tap and blot) or aspirate. Repeat two (2) additional times for a total of three (3) washes. **An automatic or manual plate washer can be used. Follow the manufacturer's instruction for proper usage. If a squeeze bottle is employed, fill each well by depressing the wash container (avoiding air bubbles) to dispense the wash. Decant the wash and repeat two (2) additional times.**
- 8. Add 0.100ml (100µl) of working substrate solution to all wells (see Reagent Preparation Section). **Always add reagents in the same order to minimize reaction time differences between wells.**  
**DO NOT SHAKE THE PLATE AFTER SUBSTRATE ADDITION**
- 9. Incubate at room temperature for fifteen (15) minutes.
- 10. Add 0.050ml (50µl) of stop solution to each well and gently mix for 15-20 seconds. **Always add reagents in the same order to minimize reaction time differences between wells.**
- 11. Read the absorbance in each well at 450nm (using a reference wavelength of 620-630nm to minimize well imperfections) in a microplate reader. **The results should be read within thirty (30) minutes of adding the stop solution.**

**10.0 CALCULATION OF RESULTS**

**A dose response curve is used to ascertain the concentration of Free Testosterone in unknown specimens within the analytical measuring range of 0.11-60 pg/ml.**

- 1. Record the absorbance obtained from the printout of the microplate reader as outlined in Example 1.
- 2. Plot the absorbance for each duplicate serum reference versus the corresponding Free Testosterone concentration in pg/ml on linear graph paper.
- 3. Connect the points with a best-fit curve.
- 4. To determine the concentration of Free Testosterone for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in pg/ml) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated). In the following example, the average absorbance 0.989 intersects the dose response curve at 2.87pg/ml Free Testosterone concentration.

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**12.1 Assay Performance**

- It is important that the time of reaction in each well is held constant to achieve reproducible results.
- Pipetting of samples should not extend beyond ten (10) minutes to avoid assay drift.
- Highly lipemic, hemolyzed or grossly contaminated specimen(s) should not be used.
- If more than one (1) plate is used, it is recommended to repeat the dose response curve.
- The addition of substrate solution initiates a kinetic reaction, which is terminated by the addition of the stop solution. The substrate and stop solution should be added in the same sequence to eliminate any time-deviation during reaction.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Failure to remove adhering solution adequately in the aspiration or decantation wash step(s) may result in poor replication and spurious results.
- Use components from the same lot. No intermixing of reagents from different batches.
- Accurate and precise pipetting, as well as following the exact time and temperature requirements prescribed are essential. Any deviation from IFU may yield inaccurate results.
- All applicable national standards, regulations and laws, including, but not limited to, good laboratory procedures, must be followed to ensure compliance and proper device usage.
- It is important to calibrate all the equipment e.g. Pipettes, Readers, Washers and/or the automated instruments used with this device, and to perform routine preventative maintenance.
- Risk Analysis- as required by CE Mark IVD Directive 98/79/EC - for this and other devices, made by Monobind, can be requested via email from [Monobind@monobind.com](mailto:Monobind@monobind.com).

**12.2 Interpretation**

- Measurements and interpretation of results must be performed by a skilled individual or trained professional.**
- Laboratory results alone are only one aspect for determining patient care and should not be the sole basis for therapy, particularly if the results conflict with other determinants.
- The reagents for the test system procedure have been formulated to eliminate maximal interference; however, potential interaction between rare serum specimens and test reagents can cause erroneous results. Heterophilic antibodies often cause these interactions and have been known to be problems for all kinds of immunoassays. (Boscato LM Stuart MC. 'Heterophilic antibodies: a problem for all immunoassays' Clin.Chem. 1988:3427-33). For diagnostic purposes, the results from this assay should be used in combination with clinical examination, patient history, and other clinical findings.
- For valid test results, adequate controls and other parameters must be within the listed ranges and assay requirements.
- If test kits are altered, such as by mixing parts of different kits, which could produce false test results, or if results are incorrectly interpreted, **Monobind shall have no liability.**
- If computer controlled data reduction is used to interpret the results of the test, it is imperative that the predicted values for the calibrators fall within 10% of the assigned concentrations.

**13.0 EXPECTED RANGES OF VALUES**

In agreement with established reference intervals for a "normal" adult population, the expected ranges for the Free Testosterone AccuBind® ELISA Test System are detailed in Table 1.

**TABLE 1**

Population	Range (in pg/ml)
Male / 20-39	9.2-34.6
Male / 40-59	6.1-30.3
Male / $\geq 60$	6.1-27.9
Female / 20-39	0.2-6.1
Female / 40-59	0.3-4.4
Female / $\geq 60$	0.5-3.4

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**14.0 PERFORMANCE CHARACTERISTICS**

**14.1 Accuracy**

The Free Testosterone AccuBind® ELISA Test System was compared with a reference ELISA method. Biological specimens from low, normal, and elevated concentrations were assayed. The total number of such specimens was 137. The least square regression equation and the correlation coefficient were computed.

**TABLE 2**

	Least Square Regression Analysis	Correlation Coefficient
Monobind (y)	$y=1.017x-0.24$	0.997
Reference (x)		

**14.2 Precision**

This study was conducted during 20 days of testing. The human serum and control sample were tested in duplicate, two times per day for a total of 80 measurements per sample. Three (3) different reagent lots, three (3) serum pools, and three (3) controls were used for the study (low, medium, and high concentration). The results of a representative lot are shown below:

**TABLE 3**

Lot	Mean (pg/ml)	Within-Run		Total	
		SD	CV%	SD	CV%
N=32					
Ctrl 1	2.51	0.09	3.7%	0.20	7.8%
Ctrl 2	10.98	0.40	3.6%	0.96	8.7%
Ctrl 3	22.72	0.83	3.6%	2.18	9.6%
Serum 1	0.98	0.06	5.9%	0.12	12.4%
Serum 2	4.53	0.26	5.7%	0.36	8.0%
Serum 3	53.62	4.24	7.9%	4.32	8.1%

**TABLE 4**

	Mean (pg/ml)	Within-Run Precision		Within-Kit Precision		Total Precision (n=80)	
		SD	CV%	SD	CV%	SD	CV%
Ctrl 1	2.48	0.11	4.57	0.20	8.20	0.21	8.51
Ctrl 2	11.04	0.47	4.23	0.84	2.60	0.87	8.00
Ctrl 3	23.24	1.00	4.31	1.80	7.73	1.83	7.95
Pnt 1	0.97	0.05	4.88	0.09	9.14	0.08	9.43
Pnt 2	4.62	0.23	4.88	0.32	6.89	0.33	7.20
Pnt 3	54.66	3.25	5.95	3.92	7.17	4.13	7.55

**14.3 Detection Limits**

The LOB (limit of the blank), the LOD (limit of detection) and the LoQ (limit of quantitation) were determined in accordance with CLSI EP 17-A A guideline, Protocols for Determination of Limits of Detection).

**TABLE 5**

LoB	LoD	LoQ
0.0295 pg/ml	0.0519 pg/ml	0.0519 pg/ml

**14.4 Cross Reactivity**

Cross reactivity was determined by testing those compounds most likely to interfere with the Monobind Free Testosterone ELISA Test System. The specificity of the assay was determined in accordance with CLSI EP07-A2. The results of the cross-reactivity study are as follows.

**TABLE 6**

Sample	Conc. (ng/ml)	Cross Reactivity	
		Spiked Serum	Blank Serum
11-Deoxycortisol	1000	0.000%	ND
11-KetoTestosterone	10	0.647%	0.519%
11 $\beta$ -Hydroxytestosterone	100	0.065%	0.054%
17 $\alpha$ -ethynyl estradiol	1000	0.000%	ND
17 $\alpha$ -Estradiol	1000	0.000%	0.000%
17 $\beta$ -Estradiol	100	0.000%	ND
17-Hydroxypregnenolone	1000	0.000%	ND
17-Hydroxprogesterone	10	0.000%	0.000%
3-EstriolGluc	1000	0.000%	ND
3-EstriolSul	1000	0.000%	ND
3 $\beta$ -Androstanediol	500	0.000%	ND
5 $\alpha$ -Dihydrotestosterone	100	0.054%	0.042%
Aldosterone	8000	0.000%	0.000%

Amitriptyl HCl	1000	0.000%	ND
Androsterone	1000	0.000%	ND
Androstenedione	1000	0.004%	0.002%
Clomiphene Citrate	1000	0.000%	ND
Corsticosterone	1000	0.000%	0.000%
Corstisone	1000	0.000%	0.000%
Cortisol	1000	0.000%	0.000%
Cyproterone acetate	1000	0.000%	ND
D-5-Androstene-3 $\beta$ ,17 $\beta$ -diol	1000	0.000%	ND
Danazol	1000	0.000%	ND
DHEA	100000	0.000%	0.000%
DHEA-S	1000	0.000%	0.000%
Desogestrel	100	0.000%	ND
Dexamethasone	1000	0.000%	ND
Epitestosterone	1000	0.001%	0.001%
Estrilol	1000	0.000%	0.000%
Estrone	1000	0.000%	0.000%
Ethisterone	1000	0.000%	0.000%
Ethinodiol	1000	0.000%	0.000%
Ethinodiol diacetate	50	0.000%	ND
Flunisolide	1000	0.000%	ND
Fluoxymesterone	1000	0.000%	ND
Lynestrol	1000	0.000%	ND
Medoxyprogesterone acetate	1000	0.000%	ND
Methyl Testosterone	100	0.000%	ND
Mestranol	1000	0.000%	ND
Norethindrone	50	0.000%	ND
Norethindrone acetate	50	0.000%	ND
Norgestimate	1000	0.000%	ND
Norgestrel (Levonorgestrel)	50	0.000%	ND
Norethynodrel	50	0.000%	ND
Oxymetholone	100	0.000%	ND
Prednisolone	1000	0.000%	ND
Prednisone	800	0.000%	0.000%
Progesterone	1000	0.000%	0.000%
Salbutamol	1000	0.000%	ND
Spirolactone	1000	0.000%	0.000%
Stanozolol	1000	0.000%	0.000%
Testosterone Cypionate	12	0.002%	0.000%
Testosterone enanthate	100	0.000%	0.000%
Testosterone SO4	1000	0.004%	0.003%
Testosterone Propionate	1000	0.000%	0.000%
Testosterone Undecanoate	12	0.011%	0.053%
Triamcinolone	50	0.000%	0.000%

**14.5 Interference**

Using CLSI-A2 Interference Testing in Clinical Chemistry as a guide, potential interferents were tested utilizing charcoal-stripped human serum spiked with known concentrations of interferent. The following results of % binding values even at higher than normal interferent levels indicate that there is no significant binding on the free testosterone-HRP conjugate.

**TABLE 7**

Substance	Highest concentration at which no significant interference was observed
Acetaminophen	20 mg/dl
Acetylcysteine	150 mg/dl
Ascorbic Acid	6 mg/dl
Bilirubin Conjugated	15 mg/dl
Bilirubin Unconjugated	20 mg/dl
Biotin	100 ng/ml
Caffeine	6 mg/dl
Cholesterol	503 mg/dl
Creatine	30 mg/dl
Dextran	5000 mg/dl
Digoxin	6.1 ng/ml
Doxycycline	50 mg/L
Erythromycin	6 mg/dl
Gentamicin	1 mg/dl
HAMA	440 ng/ml
Heparin	3 U/ml
Hemoglobin	500 mg/dl
Human Serum Albumin	2.5 g/dl
Ibuprofen	50 mg/dl
Immunoglobulin G	4 g/dl
Levodopa	20 mg/L
Lidocaine	1.2 mg/dl
Lipemia (glycerides)	1000 mg/dl
Methyldopa	20 mg/L
Nicotine	0.1 mg/dl
Phenobarbital	15 mg/dl
Protein: Total	10.5 g/dl
Rheumatoid Factor	1110 IU/ml
Salicylic Acid	60 mg/dl
SHBG	200 $\mu$ g/ml
Triglycerides	900 mg/dl
Urea	500 mg/dl

**15.0 REFERENCES**

- Dorfman, RI and Shipley, RA, ED: Androgens, New York.; John Wiley and Sons, 1956.
- Horton R, Tait JF: Androstenedione production and interconversion rates measured in peripheral blood and studies on the possible site of conversion to testosterone. J.Clin Invest 45: 301-303, 1966.

3. Faiman C and Winter. JSD, Reyes, FI, *Clin Obstet Gynaecol*, 3, 467 (1976).
4. Sizonenka, PC, *Pediatrician*, 14, 191 (1987).
5. Cummings DC, Wall SR: Non sex hormone binding globulin bound testosterone as a marker for hyperandrogenism. *J. Clin Endocrinol Metab*. 61:873-876, 1985.
6. Lashansky, G, et. al., *J Clin Endocrinol Metab*, 58, 674 (1991)
7. Tietz, NW, ED: *Clinical Guide to Laboratory Tests*, 3<sup>rd</sup> ed. Philadelphia, WA Saunders Co, 1995.
8. Smith SW: Free Testosterone, *AACC Endo*; 11(3), 59-62 (1993).

Revision: 5      Date: 2018-Oct-01      DCO: 1315  
 MP5325      Product Code: 5325-300

Size	96(A)	192(B)	
Reagent (fill)	A	1ml set	1ml set
	B	1ml set	1ml set
	C	1 (13ml)	2 (13ml)
	D	1 plate	2 plates
	E	1 (20ml)	1 (20ml)
	F	1 (7ml)	2 (7ml)
	G	1 (7ml)	2 (7ml)
	H	1 (8ml)	2 (8ml)

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### Glossary of Symbols (EN 980/ISO 15223)

In Vitro - Diagnostic Medical Device

Temperature Limitation Storage Condition (2-8° C)

Consult Instructions for Use

Catalogue Number

Contains Sufficient Test for Z

Batch Code

Used By (Expiration Day)

Date of Manufacturer

Manufacturer

Authorized Rep in European Country

European Conformity

**PERFORMANCE CHARACTERISTICS – DETECTION LIMIT STUDY**

PURPOSE	Define the low end of detection of the Free Testosterone AccuBind® ELISA Test System.
SCOPE	<p>Determining the limit of the blank (LoB), limit of detection (LoD), and the limit of quantitation (LoQ) for the Free Testosterone ELISA in conformance with CLSI EP17-A2 Protocols for Determination of Limits of Detection of Limits of Quantitation (2012).</p> <p>The limit of blank (LoB) is defined as the highest value expected to be seen in a series of observed test results on a sample that contains no analyte.</p> <p>The limit of detection (LoD) is defined as the actual concentration at which an observed test result is likely to exceed the limit of blank and may therefore be declared as a detectable value.</p> <p>The limit of quantitation (LoQ) is defined as the lowest actual concentration which the analyte is reliably detected to be acceptable for clinical use.</p>
RESPONSIBILITY	<p>The Research and Development Officer and Quality Control technician have the responsibility to prepare samples, monitor laboratory site, record and report testing data per this procedure.</p> <p>The project coordinator is assigned by the R&amp;D Manager to coordinate, monitor, and record data during the duration of the study.</p>
STUDY SITE	Tests are performed at Monobind, Inc., by personnel trained and familiar with performing the Free Testosterone ELISA test System

**MATERIALS**

1. Monobind Free Testosterone AccuBind® ELISA Test System
2. Low or non-detectable testosterone serum samples
3. Low serum testosterone samples at or below the minimum assay measuring range

**TEST PROCEDURE****A. Test Sample Preparation**

1. Process samples according to accepted laboratory practices and procedures for handling, storing, and processing biological samples.
2. Obtain sufficient Free Testosterone charcoal-stripped human serum samples.
3. Screen the sample to determine its Free Testosterone value utilizing the Free Testosterone AccuBind® ELISA Test System.

**B. Detection Studies****1. Limit of Blank (LoB)**

- a. Obtain a Free Testosterone stripped sample(s) that is considered the “lowest or analyte free” sample.
- b. Utilizing 3 samples over a 3 day period, 24 per run, 2 runs a day, for a total of 6 assays run and 144 measurements.
- c. Calculate the LoB per CLSI, EP17-A, Protocols for Determination of Limits of Detection and Limits of Quantitation for the full 144 measurements over the three samples.

$$\text{LoB} = \text{Mean}_{\text{BLANK}} + (1.645/(1-(1/(4f))))\text{SD}_{\text{BLANK}}$$

where  $f = N_s - K$ , and  $N_s =$  The number of low level sample results, and  $K =$  The number of low level samples.

**2. Limit of Detection (LoD)**

- a. Spike testosterone free serum with high grade testosterone to obtain desired high and low concentrations.
- b. Create ten (10) low level Free Testosterone samples that range from zero to four (4) times the LoB value according to the dilution ratios shown below:

**Assigned Value for low level specimens**

	1st Step Dilution		Assigned Value
	Serum (pg/L)	CalcVal after X times dil (pg/L)	Calculated Value (pg/L)
Serum-1	0.771	0.9000	0.0200
Serum-2	0.812	0.8000	0.0400
Serum-3	0.926	0.7000	0.0600
Serum-4	1.022	0.6000	0.0800
Serum-5	1.145	0.5000	0.1000
Serum-6	0.771	0.4000	0.1200
Serum-7	0.812	0.3000	0.1400
Serum-8	0.926	0.2000	0.1600
Serum-9	1.022	0.1000	0.1800
Serum-10	1.145	0.0000	0.2000

- c. Generate 4 replicates for 400 total measurements over a ten (10) day period.
- d. Calculate the LoD per the guidelines as described above in 1.d. Do not average the replicates when calculating the LoD.

$$LoD = LoB + (1.645/(1-(1/(4f))))SD_{Detection} \times 1.645$$

3. Limit of Quantitation (LoQ)

- a. If total error (TE) is < than goal for total error, then LoD = LoQ.
- b. If total error fails, assay forty (40) samples with a Free Testosterone value greater than the LoD value with the Free Testosterone ELIS Test System.

C. Data Collection and Statistical Analysis

- 1. Document the observed results.
- 2. Utilize the Analyse-It software to perform the statistical analysis.

ACCEPTANCE CRITERIA

- 1. The quality control values shall fall within the established mean and ± 2 SD for the assay to be accepted.
- 2. Total Error ≤ 0.05 pg/ml
- 3. LoB ≤ 0.1 pg/ml
- 4. LoD ≤ 0.3 pg/ml
- 5. LoQ ≥ LoD

SUMMARY

Total Error	Limit of the Blank (LoB)	Limit of Detection (LoD)	Limit of Quantitation (LoQ)
0.027 pg/ml	0.0295 pg/ml	0.0519 pg/ml	0.0519 pg/ml

## CONCLUSION

Based on the data shown, the highest value expected to be seen in a series of observed test results on a sample that contains no analyte is 0.0295 pg/ml. The limit of detection (LoD) or the actual concentration at which an observed test result is likely to exceed the limit of blank and may therefore be declared as a detectable value is 0.0519 pg/ml.

**PERFORMANCE CHARACTERISTICS – CROSS REACTIVITY & INTERFERENCE****A. Cross Reactivity**

**PURPOSE** Evaluate whether various analytes or substances cross-react with the quantitation of Free Testosterone when utilizing the Free Testosterone AccuBind® ELISA Test System

**SCOPE** Determine the cross-reactivity of various analytes or substances with the Free Testosterone ELISA test system according to CLSI EP07-A2, Interference Testing in Clinical Chemistry (2005).

**RESPONSIBILITY** The Research and Development Officer and Quality Control technician have the responsibility to prepare samples, monitor laboratory site, record and report testing data per this procedure.

The project coordinator is assigned by the R&D Manager to coordinate, monitor, and record data during the duration of the study.

**STUDY SITE** Tests are performed at Monobind, Inc., by personnel trained and familiar with performing the Free Testosterone AccuBind® ELISA test System

**MATERIALS**

1. Monobind Free Testosterone AccuBind® ELISA Test System
2. Charcoal stripped human serum sample with no free testosterone
3. Male serum with endogenous free testosterone
4. Purified cross-reactants from commercial sources (e.g. Sigma-Aldrich)

**TEST PROCEDURE****A. Sample Preparation**

1. Process samples according to accepted laboratory practices and procedures for safe handling, storing, and processing biological samples.
2. The male serum is a pool of five individual samples. Assay to determine the free testosterone value, and record.
3. The blank serum is human serum that has been charcoal-stripped of free testosterone.

**B. Test procedure**

1. The interferent stocks solutions were prepared by dissolving the interfering compounds in appropriate organic buffer.
2. Spike the male serum and blank serum samples with known concentrations of cross-reactants to commercially accepted levels as specified in Table 1.

**Table 1**

<b>Sample</b>	<b>Cross Reactant Conc. (ng/ml)</b>
11-Deoxycortisol	1000
11-KetoTestosterone	10
11 $\beta$ -Hydroxytestosterone	100
17 $\alpha$ -ethynyl estradiol	1000
17 $\alpha$ -Estradiol	1000
17 $\beta$ -Estradiol	100
17-Hydroxypregnenolone	1000
17-Hydroxprogesterone	10
3-EstriolGluc	1000
3-EstriolSul	1000
3 $\beta$ -Androstanediol	500
5 $\alpha$ -Dihydrotestosterone	100
Aldosterone	8000

Amitriptyl HCl	1000
Androsterone	1000
Androstenedione	1000
Clomiphene Citrate	1000
Corsticosterone	1000
Corstisone	1000
Cortisol	1000
Cyproterone acetate	1000
D-5-Androstene-3 $\beta$ ,17 $\beta$ -diol	1000
Danazol	1000
Dehydroepiandrosterone	100000
Dehydroepiandrosterone Sulfate	1000
Desogestrel	100
Dexamethasone	1000
Epistestosterone	1000
Estriol	1000
Estrone	1000
Ethisterone	1000
Ethinodiol	1000
Ethinodiol diacetate	50
Flunisolide	1000
Fluoxymesterone	1000
Lynestrol	1000
Medoxyprogesterone acetate	1000
Methyl Testosterone	100
Mestranol	1000
Norethindrone	50
Norethinodrone acetate	50
Norgestimate	1000
Norgestrel (Levonorgestrel)	50
Norethynodrel	50
Oxymetholone	100
Prednisolone	1000
Prednisone	800
Progesterone	1000
Salbutamol	1000
Spirolactone	1000
Stanozolol	1000
Testosterone Cyprionate	12
Testosterone Undecanoate	12
Testosterone enanthate	100
Testosterone SO <sub>4</sub>	1000
Testosterone Propionate	1000
Triamcinolone	50

- Assay the samples in duplicate with the Free Testosterone AccuBind® ELISA Test System.
- Using the average of the duplicates, calculate the cross-reactivity, ensuring the units are consistent:

$$\left( \frac{\text{Observed Value} - \text{Unspiked Value}}{\text{Mass of cross-reactant}} \right) \times 100$$

## ACCEPTANCE CRITERIA

1. The % variation from specimens tested without interfering substance should be  $\leq 10\%$
2. The quality control values shall fall within the established mean and  $\pm 2SD$  for the assay to be accepted.

## SUMMARY

Sample	Conc. (ng/ml)	Male Serum Spiked Samples		Blank Serum Spiked Samples	
		Conc (pg/ml)	Cross Reactivity	Conc. (pg/ml)	Cross Reactivity
Base Pool	0	7.408	-	0	
11-Deoxycortisol	1000	7.269	0.000%	0.000	ND
11-KetoTestosterone	10	72.139	0.647%	51.855	0.519%
11 $\beta$ -Hydroxytestosterone	100	72.815	0.065%	54.000	0.054%
17 $\alpha$ -ethynyl estradiol	1000	7.351	0.000%	0.000	ND
17 $\alpha$ -Estradiol	1000	7.282	0.000%	0.153	0.000%
17 $\beta$ -Estradiol	100	7.363	0.000%	0.000	ND
17-Hydroxypregnenolone	1000	7.35	0.000%	0.000	ND
17-Hydroxprogesterone	10	7.969	0.000%	1.012	0.000%
3-EstriolGluc	1000	7.825	0.000%	0.000	ND
3-EstriolSul	1000	11.729	0.000%	0.000	ND
3 $\beta$ -Androstanediol	500	7.917	0.000%	0.000	ND
5 $\alpha$ -Dihydrotestosterone	100	60.975	0.054%	42.382	0.042%
Aldosterone	8000	12.80	0.000%	4.919	0.000%
Amitriptyl HCl	1000	7.301	0.000%	0.000	ND
Androsterone	1000	8.061	0.000%	0.000	ND
Andronstenedione	1000	50.552	0.004%	19.248	0.002%
Clomiphene Citrate	1000	7.263	0.000%	0.000	ND
Corsticosterone	1000	10.099	0.000%	0.777	0.000%
Corstisone	1000	9.90	0.000%	0.478	0.000%
Cortisol	1000	7.316	0.000%	0.298	0.000%
Cyproterone acetate	1000	10.70	0.000%	0.000	ND
D-5-Androstene-3 $\beta$ ,17 $\beta$ -diol	1000	7.782	0.000%	0.000	ND
Danazol	1000	12.30	0.000%	0.000	ND
Dehydroepiandrosterone	100000	7.311	0.000%	7.457	0.000%
Dehydroepiandrosterone Sulfate	1000	7.45	0.000%	0.238	0.000%
Desogestrel	100	7.436	0.000%	0.000	ND
Dexamethasone	1000	7.381	0.000%	0.000	ND
Epistestosterone	1000	21.612	0.001%	8.215	0.001%
Estriol	1000	7.368	0.000%	0.122	0.000%
Estrone	1000	7.679	0.000%	0.036	0.000%
Ethisterone	1000	8.597	0.000%	0.003	0.000%
Ethinodiol	1000	8.042	0.000%	0.114	0.000%
Ethinodiol diacetate	50	7.563	0.000%	0.000	ND

Flunisolide	1000	7.456	0.000%	0.000	ND
Fluoxymesterone	1000	7.450	0.000%	0.000	ND
Lynestrol	1000	7.395	0.000%	0.000	ND
Medoxyprogesterone acetate	1000	7.426	0.000%	0.000	ND
Methyl Testosterone	100	7.163	0.000%	0.000	ND
Mestranol	1000	7.338	0.000%	0.000	ND
Norethindrone	50	7.541	0.000%	0.000	ND
Norethinodrone acetate	50	7.428	0.000%	0.000	ND
Norgestimate	1000	7.478	0.000%	0.000	ND
Norgestrel (Levonorgestrel)	50	7.463	0.000%	0.000	ND
Norethynodrel	50	7.544	0.000%	0.000	ND
Oxymetholone	100	7.401	0.000%	0.000	ND
Prednisolone	1000	7.828	0.000%	0.000	ND
Prednisone	800	10.442	0.000%	0.966	0.000%
Progesterone	1000	11.025	0.000%	0.248	0.000%
Salbutamol	1000	8.026	0.000%	0.000	ND
Spironolactone	1000	7.310	0.000%	0.112	0.000%
Stanozolol	1000	7.442	0.000%	0.081	0.000%
Testosterone enanthate	100	7.304	0.000%	0.044	0.000%
Testosterone SO4	1000	75.970	0.004%	28.182	0.003%
Testosterone Propionate	1000	9.589	0.000%	1.008	0.000%
Triamcinolone	50	7.544	0.000%	0.043	0.000%
Base Pool for below analytes		38.436	-	0.077	-
Testosterone Cypionate*	12	38.685	0.002%	0.021	0.000%
Testosterone Undecanoate*	12	38.262	-0.001%	0.015	-0.001%

\*C<sub>MAX</sub> for the two compounds is 11 ng/ml or 11000 pg/ml

#### CONCLUSION

1. No cross reactivity above 0.647% was detected with the performance of Free Testosterone AccuBind® ELISA Test System upon addition of massive amounts of selected substances to the human serum pool.

**B. Interference**

**PURPOSE** Conduct interference testing to identify potential interferents when utilizing the Free Testosterone ELISA Test System

**SCOPE** Determine potential intereferents with the Free Testosterone ELISA test system according to CLSI EP07-A2, Interference Testing in Clinical Chemistry (2005)

**RESPONSIBILITY** The Research and Development Officer and Quality Control technician have the responsibility to prepare samples, monitor laboratory site, record and report testing data per this procedure.

The project coordinator is assigned by the R&D Manager to coordinate, monitor, and record data during the duration of the study.

**STUDY SITE** Tests are performed at Monobind, Inc., by personnel trained and familiar with performing the Free Testosterone AccuBind® ELISA test System

**MATERIALS**

1. Monobind Free Testosterone AccuBind® ELISA Test System
2. Charcoal-stripped human serum pool
3. Commercially available interferents for spiking

**TEST PROCEDURE****A. Sample Preparation**

1. Process samples according to accepted laboratory practices and procedures for safe handling, storing, and processing biological samples.
2. The male serum is a pool of five individual samples. Assay to determine the free testosterone value, and record.
3. The blank serum is human serum that has been charcoal-stripped of free testosterone.

**B. Test Procedure**

1. Create a base pool containing the sample matrix of interest, which can comprise a volume of 4-6 ml.
2. The interferent stocks solutions were prepared by dissolving the interfering compounds in appropriate buffer.
3. It is recommended that screening experiments be done with both low and high analyte pools, as interference may be dependent upon analyte concentration.
4. Spike human serum samples with increasing amounts of interferent. The interferent volume should be no more than 5% of the matrix volume.
5. Document the concentrations of interferent added.
6. Assay according to instructions for Free Testosterone AccuBind® ELISA Test System.
7. Determine if interferent substances cause intereference at specified concentrations by calculating percent interference.

**C. Data Collection and Statistical Analysis**

1. Document the observed results.
2. Observe for effect of interferent on a sample by comparing the concentration of unspiked to spiked.
3. Using the average of the duplicates, calculate the percent interference, ensuring the units are consistent:

$$\frac{\text{Observed Value} - \text{Unspiked Value}}{\text{Unspiked Value}}$$

ACCEPTANCE CRITERIA

1. Non-significant interference:
  - a. The % variation from specimens tested without interfering substance should be  $\leq 10\%$
  - b. Recovery after the addition of interferents should be  $\geq 90\%$
2. The quality control values shall fall within the established mean and  $\pm 2SD$  for the assay to be accepted.

SUMMARY

Effect of Human Serum Albumin (HSA)

Sample	Concentration Added (g/dL)			
	0	1.25	2.5	5
1	2.298	1.993	2.220	1.997
2	10.832	10.606	9.254	8.357
3	5.018	4.747	4.513	4.182
4	0.989	1.088	1.087	1.086

Effect of Sex Hormone Binding Globulin (SHBG)

Sample ID	SHBG Added	OD	Percent
	$\mu\text{g/ml}$	450nm	B/B(0)
1	0	2.487	100.00%
2	6.25	2.473	99.42%
3	12.5	2.483	99.82%
4	25	2.458	98.81%
5	50	2.460	98.91%
6	100	2.452	98.59%
7	200	2.436	97.93%

		Male Serum High Spiked Samples Base Concentration: 38.498 pg/ml	Male Serum Mid-Level Spiked Samples Base Concentration: 7.916 pg/ml	Blank Serum Spiked Samples
Sample	Concentration Added	Interference	Interference	Interference
Acetaminophen	20 mg/dl	0.699%	2.804%	ND
Acetylcysteine	150 mg/dl	1.597%	-0.189%	ND
Ascorbic Acid	6 mg/dl	3.720%	1.832%	ND
Bilirubin Conjugated	15 mg/dl	2.268%	-2.830%	ND
Bilirubin Unconjugated	20 mg/dl	1.626%	-2.893%	ND
Biotin	100 ng/ml	3.704%	-2.590%	ND
Caffeine	6 mg/dl	2.151%	1.983%	ND

Cholesterol	503 mg/dl	2.551%	4.156%	ND
Creatine	30 mg/dl	1.574%	3.550%	ND
Dextran	5000 mg/dl	0.899%	2.665%	ND
Digoxin	6.1 ng/ml	-1.735%	-3.499%	ND
Doxycycline	50 mg/L	-1.995%	-4.775%	ND
Erythromycin	6 mg/dl	0.714%	-0.076%	ND
Gentamicin	1 mg/dl	-1.998%	-4.207%	ND
Hemoglobin*	500 mg/dl	0.774%	-2.700%	ND
Heparin	3 U/ml	-3.008%	-5.584%	ND
Ibuprofen	50 mg/dl	-1.961%	0.227%	ND
Immunoglobulin G	4 g/dl	-1.397%	-1.592%	ND
Levodopa	20 mg/L	3.075%	-1.187%	ND
Lidocaine	1.2 mg/dl	0.158%	1.377%	ND
Lipemia	1000 mg/dl	-0.696%	-3.209%	ND
Methyldopa	20 mg/L	-0.722%	-3.095%	ND
Nicotine	0.1 mg/dl	-6.003%	0.417%	ND
Phenobarbital	15 mg/dl	-7.374%	-0.644%	ND
Protein: Total	10.5 g/dl	2.400%	-2.994%	ND
Rheumatoid Factor	1110 IU/ml	-1.839%	-4.775%	ND
Salicylic Acid	60 mg/dl	-5.117%	4.005%	ND
SHBG	96nM	-4.592%	-2.451%	ND
Triglycerides	900 mg/dl	-3.372%	-4.169%	ND
Urea	500 mg/dl	-1.704%	0.733%	ND

\*Hemoglobin interference of male serum mid-spiked samples were tested with a base pool concentration of 9.409pg/ml

Human anti-mouse antibodies (HAMA) interference was evaluated with human serum samples with low and high endogenous HAMA. Less than 10% interference was observed for HAMA at 440 ng/ml and 185 ng/ml.

#### CONCLUSION

1. The results showed % binding values between 97.93% – 99.42% even at higher than normal interferent levels for HSA and SHBG. Results showed no interference above 10% for interferents tested on high, mid-level, and blank Free Testosterone concentration. In conclusion, the results showed that there was no significant binding by interferents on the free testosterone-HRP conjugate