



GE HealthCare

Versana Premier/Versana Premier Lotus User Manual



5928218-1EN Rev. 6

Version R3

GENERAL USER DOCUMENTATION

© 2023 By GE HealthCare



GE HealthCare

Regulatory Requirement

The Versana Premier/Versana Premier Lotus complies with regulatory requirements of the following European Regulation 2017/745 EU concerning medical devices.



First CE Marked in 2023.

This User Manual is a reference for the Versana Premier/Versana Premier Lotus. It applies to Version R3 software for the Versana Premier/Versana Premier Lotus ultrasound system.

GE Healthcare
P.O. Box 414, Milwaukee, Wisconsin 53201 U.S.A.
(Asia, Pacific, Latin America, North America)

GE Healthcare GmbH
Beethovenstrasse 239
Postfach 11 05 60
D-42655 Solingen GERMANY
TEL: 49 212.28.02.208; FAX: 49 212.28.02.380

Revision History

Table -1 Reason for Change

Revision	Date (YYYY-MM-DD)	Reason for change
Rev. 1	2023-10-19	Initial release.
Rev. 2	2023-10-30	Add India factory address.
Rev. 3	2024-02-21	• Update first CE mark year. • Remove Tricify Uplink. • Update graphics and description.
Rev. 4	2024-10-09	Update InSite ExC service content.
Rev. 5	2024-11-20	Update Digital Expert Connect section.
Rev. 6	2025-03-06	Add Full Screen and update Worklist description.

Please verify that you are using the latest revision of this document. Information pertaining to this document is maintained on MyWorkshop. If you need to know the latest revision, contact your distributor, local GE HealthCare Sales Representative or in the USA call the GE HealthCare Ultrasound Clinical Answer Center at 1 800 682 5327 or 1 262 524 5698.

Regulatory Requirements

Conformance Standards

The following classifications are in accordance with the IEC/EN 60601-1:6.8.1:

- According to European Medical Device Regulation 2017/745, this is Class IIa Medical Device.
- According to IEC 60601-1:
Equipment is Class I or Internally Powered equipment, Type BF or Type CF Applied Part for probes, DEFIBRILLATION-PROOF Type CF Applied Part for ECG.
Not suitable to be used in an OXYGEN RICH ENVIRONMENT.
Continuous Operation.
- According to CISPR 11:
Equipment is Group 1, Class A ISM Equipment.
- According to IEC 60529:
The footswitch rate IPx8 is suitable for use in surgical rooms.
Probe head (immersible portion) and cable are IPX7. Probe connector is not waterproof.

This product complies with the regulatory requirement of the following:

- Council Regulation 2017/745 concerning medical devices: the CE label affixed to the product testifies compliance to the Regulation.

The location of the CE marking is shown in the Safety chapter of this manual.



Authorized EU Representative

European registered place of business:

GE Medical Systems SCS

283 rue de la Minière

78530 BUC, France



Authorized Swiss Representative

GE Medical Systems (Schweiz) AG

Europa-Strasse 31

8152 Glattbrugg

Switzerland

- International Electrotechnical Commission (IEC).

- IEC/EN 60601-1 Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance.
- IEC/EN 60601-1-2 Medical Electrical Equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances- Requirements and tests.
- IEC/EN 60601-1-6 (Usability), EN ISO 20417 (Information supplied with medical devices).
- IEC/EN 60601-2-37 Medical electrical equipment - Part 2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic monitoring equipment.
- IEC 62359 Ultrasonics - Field Characterization - Test methods for the determination of thermal and mechanical indices related to medical diagnostic ultrasound fields.
- International Organization of Standards (ISO)
ISO 10993-1 Biological evaluation of medical devices.
- ANSI/AAMI ES60601-1 Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance.
- Canadian Standards Association (CSA).
CAN/CSA-22.2, NO. 60601-1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.
- CSA 22.2, 601.1 Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance.
- Medical Device Good Manufacturing Practice Manual issued by the FDA (Food and Drug Administration, Department of Health, USA).

Certifications

General Electric Medical Systems (China) Co., Ltd. is ISO 9001 and ISO 13485 certified.

Original documentation

- The original document was written in English.

Country-Specific Approval

Brazil

Registro ANVISA N°: 80071260409

Authorized Representative in Kazakhstan

Authorized representative in Kazakhstan which is responsible for accepting claims is GE HealthCare Kazakhstan LLP.

English	Kazakh	Russian
GE HealthCare Kazakhstan LLP 26/41, Zenkova Street, Medeu District, Almaty, 050010, Kazakhstan T +7 727 3560020	«ДжиИ Хэлскеа Қазақстан» ЖШС Қазақстан, Алматы қаласы, Медеу ауданы, көшесі ЗЕНКОВ, үй 26/41, пошталық индексі 050010 T +7 727 3560020	ТОО «ДжиИ Хэлскеа Казахстан» Казахстан, город Алматы, Медеуский район, улица Зенкова, дом 26/41, почтовый индекс 050010 T +7 727 3560020

Authorized Representative in Ukraine

Ультразвукова діагностична система Versana Premier / Versana Premier Lotus	
 UA.TR.116	Знак відповідності технічним регламентам
	GE Medical Systems (China) Co., Ltd. No.19, Changjiang Road, Wuxi National Hi-Tech Development Zone, 214028 Jiangsu, P.R. China ДЖИИ Медікал Системс (Чайна) Ко., Лтд. №. 19, Чангжіанг Родд, Вуксі Нашіонал Хай-Тек Дев. Зона 214028 Джангсу, Китай
Уповноважений представник в Україні: ТОВ «ДЖИИ ХЕЛСКЕР УКРАЇНА.» 02000, Україна, місто Київ, проспект Павла Тичини, будинок 1В, офіс 408 Тел.: +38 (044) 390 06 08	
Дата останнього перегляду посібника користувача є датою останньої редакції та вказана у посібнику користувача.	

Importer Information

Türkiye

GE Medical Systems Türkiye Ltd. Şti.
Esentepe Mah. Harman Sok. No: 8
34394 Şişli İstanbul Türkiye

Brazil

INFORMAÇÕES DO DETENTOR DA NOTIFICAÇÃO NO BRASIL
GE HEALTHCARE DO BRASIL COMÉRCIO E SERVIÇOS PARA
EQUIPAMENTOS MEDICOS-HOSPITALARES LTDA.

Av. Magalhães de Castro, 4800 – Andar 10 Conj. 101 e 102,
Torre 3 - Cidade Jardim - CEP: 05676-120 - São Paulo/SP – Brasil

CNPJ: 00.029.372/0001-40

Table of Contents

Chapter 1 - Getting Started

Overview	12
Attention	12
Documentation	13
Principles of Operation	14
Contraindication	14
Patient population	14
Clinical Benefit	15
Prescription Device	15
Intended Use	15
Preparing the System for Use	16
Console Overview	19
Peripheral/Accessory Connection	20
Battery	21
Gel Warmer (Option)	26
Wired Footswitch (Option)	26
Control Panel Map	28
Indicator Light (For Versana Premier Lotus only)	30
Physical A/N keyboard (Option)	30
Top/Sub Menu	31
Touch Panel	32
Monitor	35
Moving the system	38
Before moving the system	38
Wheels	39
Moving the System	40
Transporting the System	41
System Start-Up	42
Connecting the System	42
Power On	44

Table of Contents

Power Off	47
Log Off	47
Standby Mode	48
Crash Recovery Instructions	48
Probes	50
Connecting the Probe	50
Cable handling	51
Activating/Deactivating the Probe	51
Disconnecting the Probe	51
Transporting Probes	52
Storing the Probe	52
Preparing for an Exam	54
Patient Screen	54
Scanning a New Patient	62
Entering a Patient List	64
Starting a new exam on an existing patient	65
Scanning without entering any patient data	65
Ending a Patient	65
Changing Current Patient to Existing Patient (with Patient ID)	66
Changing Current Patient to Existing Patient (without Patient ID)	66
Delete the existing patient/exam/image	67
Barcode Reader	68

Chapter 2 - Performing an Exam

Optimizing the Image	74
Optimizing B-Mode	74
Optimizing M-Mode	80
Optimizing Color Flow Mode	83
Optimizing M Color Flow	86
Optimizing Doppler Mode	86
Other Controls	90
Split screen	93
Freezing an image	93

Activating CINE Mode	93
Body Patterns	93
Annotating an Image	94
Measurement and Analysis	95
Location of Measurement Controls	95
B-Mode Measurements	96
M-Mode Measurements	100
Doppler Mode Measurements	102
Application M&A	108
Wide Dual Screen Measurements	145
Viewing and Editing Worksheets	146
Clinical Measurement Accuracy	147

Chapter 3 - After The Exam Is Over

Preset	152
System Presets	152
Imaging Presets	153
Storing Images and Cineloops	161
Data Transfer	163
Image/Data Management	171
Configuring Connectivity	186
Printing Options	189
Electronic Documentation (Not applicable for China)	193
Documentation Distribution	193
Using Online Help Via Help	193
Electronic media	194
System Care and Maintenance	196
Expected Service Life Description	196
Maintenance Schedule	196
Cleaning and Disinfecting the system	198
Prevention of static electricity interference	205
Disposal	205
Troubleshooting	206

Table of Contents

InSite ExC	207
eDelivery - Software update	215
Assistance	225
Quality Assurance	233
Typical Tests to Perform	233
Baselines	234
Periodic Checks	235
Results	235
System Setup	235
Test Procedures	236
Setting up a Record Keeping System	241
Ultrasound Quality Assurance Checklist	242
Probe Check	242
Image Quality Check (IQC)	248
Privacy and Security	249
Potential Hazardous Situations Resulting from Failures of the IT Network . . .	249
Warning	249
Contact Information	250
Contacting GE HealthCare Ultrasound	250
Manufacturer	254
Factory Sites	255

Chapter 4 - Safety

Owner Responsibility	258
Owner requirements	258
Notice against user modification	258
Bioeffect Disclosure	259
Training	259
ALARA	259
Safety Precautions	260
Precaution Levels	260
Hazard Symbols	260
Patient Safety Related Hazards	262

Equipment and Personnel Safety	264
Electrical safety	269
EMC (Electromagnetic Compatibility)	270
Patient Environmental Devices	281
Acoustic Output	282
RoHS Hazardous Substances	284
WEEE Passport	286
Safe Product and Packaging Disposal	286
General Caution	286
Device Labels	287
Label Icon Description	287
Label Locations	292
Probe Label Explanation	294
Probe Box Label	296

Chapter 5 - Probes and Biopsy

Probe Overview	298
Ergonomics	298
Cable handling	298
Probe orientation	298
Labeling	299
E7C8L-RS Probe Composition	299
L4-20t-RS Probe Configuration	299
Probe Usage	300
Care and Maintenance	300
Covering the Transducer using a Sterile, Protective Sheath	315
E7C8L-RS Operation	316
Probe Safety	319
Special handling instructions	322
Probe handling and infection control	323
Special Reprocessing Instructions for E7C8L-RS	324
Coupling gels	327
Probe Discussion	329

Table of Contents

Application	329
Features	330
Specifications	333
Slice Thickness Specification	334
Probe Illustration	335
Biopsy Special Concerns	339
Precautions Concerning the Use of Biopsy Procedures	339
Preparing for a Biopsy	341
Displaying the Guide zone	341
Preparing the Biopsy Guide Attachment	343
4D Probe Biopsy Needle Path Selection	351
Biopsy Needle Path Verification	351
The Biopsy Procedure	352
Post Biopsy	352

Chapter 6 - Advanced Features

Whizz Easy Style	355
Activating	355
Whizz CF	356
Using Whizz CF	356
B Flow / B Flow Color (Option)	358
B Flow	358
B-Flow Color (BFC)	360
Stress Echo (Option)	361
Getting started with a stress study	361
Image acquisition	363
Continuous Capture mode	365
Post Acquisition Features	371
Analysis	372
Editing/Creating template	374
Wall Motion Segment Setup	380
Utility Application Settings for Protocol	381
Reports	381

ECG Option	383
ECG connecting	384
To activate ECG Option	384
Operation	385
Patient preparation	386
Cleaning and Disinfection	386
Troubleshooting - Noisy ECG Signal	388
ECG Trace Monitor Display	389
ECG menu	389
Using 3D (Option)	392
Easy 3D	392
3D Acquisition Parameter Description	396
Advanced 3D	400
Movie 3D	403
Elastography (Option)	404
Activating Elastography	404
How to Use	405
Additional Notes for Elastography	406
Contrast Imaging (Option)	408
Mode	409
Parameters	409
Relationship with other controls	413
Time Intensity Curve (TIC) Analysis	414
Needle Recognition (Option)	438
Adjusting	438
Breast Care	443
Prepare the patient information and scanning	443
Breast Productivity Package	446
Breast Lesion M&A	446
BI-RADS (ACR)	447
Worksheet and Summary Worksheets	450
Breast Measure Assistant (Auto Contour)	450
Thyroid Productivity Package	452
TI-RADS (ACR)	455

Table of Contents

Worksheet and Summary Worksheets	458
Measure Assistant Thyroid (Auto Contour)	459
Scan Assistant	460
Scan Assistant Manager	460
Activating the Scan Assistant Creator	466
Exporting Scan Assistant Programs	466
My Trainer	468
Follow-up Tool	470
Follow-up tool configuration	470
Follow-up Tool Workflow	472
Worksheet	474
Report	474
QAnalysis (Option)	476
Activating QAnalysis	476
QAnalysis Screen Description	477
Selecting QAnalysis Image Range	478
Generating a Trace	479
Manual tracking of the sample area (dynamic anchored sample area)	479
Delete a trace	480
Manipulating the Sample Area	480
Disabling/Enabling the frame	483
Smoothing	484
Horizontal Sweep	484
Drift Compensation	484
Statistics	485
Trace Measurements	485
Exporting Traces (Saving the Trace Data)	485
Annotating the QAnalysis Data	487
Exiting QAnalysis	487
Whizz Label (Not applicable for China)	488
Using Whizz Label	489
Advanced urology procedures	492
ERB7/ERB Biopsy Attachment for E7C8L-RS	492
Biopsy Guideline Type Preset Selection	493

E7C8L-RS Grid Biopsy Accessories	495
Calibration Procedure	497
Stepper Volume Measurement	500
Strain Imaging	503
Using Strain	505
Optimizing Strain	505
Strain rate	506
Using Strain rate	507
Optimizing Strain rate	507
Using 4D	508
Features supported with 4D	508
4D Principles of Operation	508
4D Operational Controls	515
4D Presets	516
Performing a 4D Scan	517
Tomographic Ultrasound Imaging (TUI)	531
Whizz RenderLive	533
Whizz Follicle	535
Using Whizz Follicle	537
Whizz Follicle Report	538
Whizz Note	540
VCO for customer	544
Outside Hospital version	544
Inside Hospital version	553
Using VCO for Customer - EClient Side	557
Digital Expert Connect	558
Intensity Ratio	559
Purpose	559
Probes and presets	559
Scan Mode	559
Activation	559
IR Calculation	559
IR Calculation with Follow Up Tool	562
Worksheet and Report	563

Table of Contents

US LI-RADS (ACR)	564
US Category	564
US Visualization Score	564
User Requirement	565
Whizz Report	568
Creating a report	568
Activating the Report	568
Selecting another template	569
Factory Templates	569
Whizz Report Editor	570
Import Report Template to Ultrasound System	580
Scan Coach	582
Scan Coach Manager	582
Activating the Scan Coach Creator	588
Scan Coach Creator	589

Chapter 1

Getting Started

In this section

Overview.	12
Moving the system.	38
System Start-Up.	42
Probes.	50
Preparing for an Exam.	54

Overview

Attention

This manual contains information to operate the system safely. Advanced equipment training may be provided by a factory trained Applications Specialist for the agreed-upon time period.

Read and understand all instructions in this manual before attempting to use the ultrasound system.

Keep this manual with the equipment at all times. Periodically review the procedures for operation and safety precautions.

NOTE

The Online Help offers a quick way for the user to access the manual.

Disregarding information on safety is considered abnormal use.

Not all features, products, probes or peripherals described in this document may be available or cleared for sale in all markets. Please contact your local GE HealthCare Ultrasound representative to get the latest information.

Table 1-1 Product model configuration

Model	3P P	4P P	5P P	Fixed arm	Flexible arm	21.5 inch LCD	23.8 inch LCD	Indicator Light	Strain	Strain Rate
Versana Premier VS	✓	O	X	✓	O	✓	X	X	X	X
Versana Premier VA	X	✓	X	X	✓	X	✓	X	X	X
Versana Premier Lotus 4PP	X	✓	X	X	✓	X	✓	✓	O	O
Versana Premier Lotus 5PP	X	X	✓	X	✓	X	✓	✓	O	O
NOTE ✓: The feature/configuration is available for current model. X: The feature/configuration is not available for current model. O: The feature/configuration is optional for current model.										

NOTE

Please note that orders are based on the individually agreed upon specifications and may not contain all features listed in this manual.

NOTE

All references to standards / regulations and their revisions are valid for the time of publication of the user manual.

Documentation

**CAUTION**

Safety instructions must be reviewed before operating the unit.

Versana Premier/Versana Premier Lotus documentation is provided electronically on the ultrasound system by pressing **Help**, by viewing the documentation media via the Service Desktop, or may be ordered on paper (via H-Cat). Additionally, an electronic Instructions For Use Kit (eIFU Kit) is shipped with the ultrasound system. The eIFU Kit provides all of the manuals in all languages.

Table 1-2 Documentation

Publication	Translated	Available via Help	Available via Media in the eIFU Kit	Available on Paper (if purchased)
User Manual Provides information needed by the user to operate the system safely.	Yes (all required languages)	Yes	Yes	Yes
Advanced Reference Manual Provides Acoustic Output Data and System Measurement and Analysis Tables and Formulas.	English and French	Yes	Yes	Yes
Basic Service Manual Supplies block diagrams, spare parts, adjustments, instructions to help qualified technical personnel repair the system.	No	No	Yes	Yes
Release Notes Provides precautions and instructions that supplement the User Manual.	Yes	No	Yes	Yes

NOTE

An AIUM Booklet is shipped with systems shipped in the United States and Canada.

NOTE

The screen graphics in this manual are only for illustrational purposes. Actual screen output may differ.

Principles of Operation

Medical ultrasound images are created by computer and digital memory from the transmission and reception of mechanical high-frequency waves applied through a transducer. The mechanical ultrasound waves spread through the body, producing an echo where density changes occur. For example, in the case of human tissue, an echo is created where a signal passes from an adipose tissue (fat) region to a muscular tissue region. The echoes return to the transducer where they are converted back into electrical signals.

These echo signals are highly amplified and processed by several analog and digital circuits having filters with many frequency and time response options, transforming the high-frequency electrical signals into a series of digital image signals which are stored in memory. Once in memory, the image can be displayed in real-time on the image monitor. All signal transmission, reception and processing characteristics are controlled by the main computer. By selection from the system control panel, the user can alter the characteristics and features of the system, allowing a wide range of uses, from obstetrics to peripheral vascular examinations.

Transducers are accurate, solid-state devices, providing multiple image formats. The digital design and use of solid-state components provides highly stable and consistent imaging performance with minimal required maintenance. Sophisticated design with computer control offers a system with extensive features and functions which is user-friendly and easy to use.

Contraindication

The ultrasound system is not intended for ophthalmic use or any use causing the acoustic beam to pass through the eye.

Patient population

- Age: all ages (including embryos and fetuses)
- Location: worldwide
- Sex: male and female
- Weight: all weight categories

NOTE

Extreme obesity may change the diagnostic information ability of the device.

- Height: no limitations

Clinical Benefit

The clinical benefit of a diagnostic ultrasound device is to help healthcare professionals provide accurate diagnostic information (visualize human tissue/internal structure) that enhances the diagnostic and treatment care pathways of the patient for a variety of diseases and conditions.



CAUTION

This machine should be used in compliance with law. Some jurisdictions restrict certain uses, such as gender determination.

Prescription Device



CAUTION: United States law restricts this device to sale or use by, or on the order of a physician.

Intended Use

The Versana Premier/Versana Premier Lotus is a general-purpose diagnostic ultrasound system intended for use by qualified and trained healthcare professionals for ultrasound imaging, measurement, display and analysis of the human body and fluid.

Indications for Use

The Versana Premier/Versana Premier Lotus is a general-purpose diagnostic ultrasound system intended for use by qualified and trained healthcare professionals for ultrasound imaging, measurement, display and analysis of the human body and fluid. Versana Premier/Versana Premier Lotus clinical applications include: Fetal/Obstetrics, Abdominal, Gynecology, Urology, Pediatric, Small Parts (includes breast, testes, thyroid), Cardiac Adult, Cardiac Pediatric, Vascular/Peripheral Vascular, Musculoskeletal Conventional, Musculoskeletal Superficial, Thoracic/Pleural, Transcranial, Transrectal, Transvaginal, Transesophageal, Interventional guidance (includes tissue biopsy, fluid drainage, vascular and non-vascular access).

Modes of operation include: B, M, PW Doppler, CW Doppler, Color Doppler, Color M Doppler, Power Doppler, Harmonic Imaging, Coded Pulse, 3D/4D Imaging mode and Combined modes: B/M, B/Color, B/PWD, B/Color/PWD, B/Power/PWD.

Versana Premier/Versana Premier Lotus is intended to be used in a hospital or medical clinic.

Frequency of Use

Daily (Typically 8 hours)

Type of use

Multiple patient multiple use.

Operator Profile

- Qualified and trained Healthcare professionals, including physicians, sonographers and equivalent/comparable professions, with at least basic ultrasound knowledge.
- The operator must have read and understood the user manual.

NOTE

Only qualified physicians or sonographers should perform ultrasound scanning on human subjects for medical diagnostic reasons. Request training, if needed.

Preparing the System for Use

Site Requirements

Do not attempt to install the unit alone. General Electric, Affiliate, or Distributor Field Engineers and Application Specialists will install and setup the system.

The ultrasound system does not contain any operator serviceable internal components. Ensure that unauthorized personnel do not tamper with the unit.

Perform regular preventive maintenance.

Maintain a clean environment. Turn off, and if possible, disconnect the system before cleaning the unit.

Never set liquids on the unit to ensure that liquid does not drip into the control panel or unit.

Before the system arrives

The ultrasound unit must operate within the proper environment and in accordance with the requirements described in this section. Before using the system, ensure that the requirements are met.

Power Requirements

- A separate power outlet with a 6.5 amp Supply Mains Switch.
- Frequency: 50/60 Hz
- 100 - 240V AC ($\pm 10\%$)

Fuse Type and Nominal Value

Type: 0215010.HXP

Factory Tag: LITTELFUSE

Nominal Value: 10A, 250VAC

NOTE

For changing the fuse, please contact GE HealthCare Service.

Electromagnetic interferences

This medical equipment is approved, in terms of the prevention of radio wave interference, to be used in hospitals, clinics and other institutions which are environmentally qualified.

The use of this equipment in an inappropriate environment may cause some electronic interference to radios and televisions around the equipment.

Ensure that the following is provided for the new system:

- Take precautions to ensure that the console is protected from electromagnetic interference.

Precautions include:

- Operate the console at least 15 feet away from motors, typewriters, elevators, and other sources of strong electromagnetic radiation.
- Operation in an enclosed area (wood, plaster or concrete walls, floors and ceilings) helps prevent electromagnetic interference.
- Special shielding may be required if the console is to be operated in the vicinity of radio broadcast equipment.



CAUTION

Do not operate the system in the vicinity of a heat source, of strong electric or magnetic fields (close to a transformer), or near instruments generating high-frequency signals, such as HF surgery. These can adversely affect the ultrasound images.

Environmental Requirements

The system should be operated, stored, or transported within the parameters outlined below. Either its operational environment must be constantly maintained or the unit must be turned off.

NOTE

You may get an overheating message with regard to fan speed. Ensure adequate system/room ventilation.

Table 1-3 System Environmental Requirements (including probes)

Element	Operational (with probe)	Storage	Transport
Temperature	3 ~ 40°C (without battery) 37.4 ~ 104°F (without battery) 10 ~ 30°C (with battery) 50 ~ 86°F (with battery)	-20 ~ 50°C (without probes) -4 ~ 122°F (without probes) -5 ~ 50°C (with probes) 23 ~ 122°F (with probes)	-20 ~ 50°C (without probes) -4 ~ 122°F (without probes) -5 ~ 50°C (with probes) 23 ~ 122°F (with probes)
Humidity	30% ~ 85% non-condensing	10% ~ 90% non-condensing	10% ~ 90% non-condensing
Pressure	700 ~ 1060hPa	700 ~ 1060hPa	700 ~ 1060hPa

NOTE

The operational temperature should be 18-30 degree while the system is connecting with RAB2-6-RS.

**CAUTION**

The system cannot be used in OXYGEN rich environment.

**CAUTION**

Ensure that the probe face temperature does not exceed the normal operation temperature range.

Operating Environment

Ensure that there is sufficient air flow around the ultrasound unit when installed in a fixed location.

**CAUTION**

Do not cover the ventilation holes of the ultrasound system.

**CAUTION**

The ultrasound system and probe connector are NOT waterproof. Do not expose the device to water or any kind of liquid.

Console Overview

NOTE

The graphics of the ultrasound system in this manual are for illustration purpose only, real system may vary with different configuration.

Figure 1-1 Ultrasound system overview - example



- | | |
|----------------------------------|--------------------------------------|
| 1 LCD Monitor | 9 Probe Connector |
| 2 Touch Panel | 10 Speaker |
| 3 Gel Bottle Holder/Gel Warmer | 11 Rear Handle |
| 4 Control Panel | 12 Rear Panel |
| 5 Probe Holder | 13 Wheel |
| 6 Flexible Cable Hook/Cable Hook | 14 Flexible Arm / Fixed Arm |
| 7 Printer & DVD Box | 15 Back side Accessory Tray (Option) |
| 8 Body Left and Right Cover | 16 Back side Accessory Box (Option) |

NOTE

The console graphics in this manual are only for illustrational purposes. Refer to the system for actual appearance.

Peripheral/Accessory Connection

The ultrasound system peripherals and accessories can be properly connected using the rear connector panel.



CAUTION

For compatibility reasons, use only GE HealthCare approved probes, peripherals or accessories. **DO NOT** connect any probes, peripherals or accessories without approval by GE HealthCare.



CAUTION

The connection of equipment or transmission networks other than as specified in the user instructions can result in electric shock hazard. Alternate connections will require verification of compatibility and conformity to IEC/EN 60601-1 by the installer.



CAUTION

DO NOT touch the patient and any of the connectors on the ultrasound unit simultaneously, including ultrasound probe connectors. DO NOT touch the conducting parts of the USB, Ethernet, Video, Audio cables when connecting equipment to the unit.



CAUTION

When using peripheral device, observe all warnings and cautions given in Peripheral manufacture's manuals.



WARNING

Peripheral devices that use their own AC power source can not be attached to system.

Figure 1-2 Peripheral/Accessory Connector Panel




Table 1-4 Peripheral/Accessory Connector Panel

No.	Electronic interface	Description (Interface standards / Data rate)
1	Diagnostic LED 1, 2, 3	System Diagnostic LED
2	Ethernet port	RJ45 network Connector, IEEE 802.3 10 Base-T, 100 Base-TX and 1000 Base-T
3	USB port	Type A USB port, USB 3.0, Max 5.0 Gbps, DC5V/500mA
4	HDMI port	Video/Audio port, HDMI, 1080P
5	Type-C USB port	Type-C USB port, USB 3.0, Max 3.0, Max 5.0 Gbps, DC5V/500mA
6	Equalization conductor	The terminal to be used for connecting equipotential conductors when interconnecting (grounding) with other equipment.
7	Supply Mains Switch	Supply Mains Switch
8	AC Inlet	AC 100-240V, IEC C13 socket

Battery

The lithium ion battery provides power when an AC power source is not available. A battery in the battery bay is optional for the ultrasound system. Lithium ion batteries last longer than conventional batteries and do not require replacement as often. You can expect minimum 20 minutes of battery life with a full charged battery in use to supply power to the system.

NOTE

While scanning with the battery supplying power only, the battery life may be shorter. Always archive the data and keep your attention on the battery status. When the battery power is low, charge by connecting to the power source immediately so that scanning is not interrupted as data could be lost due to the automatic shutdown of the system.

The lithium ion technology used in your system's battery is significantly less hazardous to the environment than the lithium metal technology used in some other batteries (such as watch batteries). Used batteries should not be placed with common household waste products. Contact local authorities for the location of a chemical waste collection program nearest you.

NOTE

The battery is designed to work with this ultrasound systems only. Only use GE HealthCare recognized batteries.

Where the integrity of external protective conductor in the installation or its arrangement is in doubt, EQUIPMENT shall be operated from its battery.



WARNING

For safety of the battery, do not expose it to temperature over 60°C (140°F). Keep it away from fire and other heat sources.

**CAUTION**

The battery may not charge properly if the environmental temperature is below 10°C (50°F) or above 35°C(95°F) which may disrupt the power supplied to the system.

Failure to read and follow the battery pack safety instructions and warnings below may result in personal injury, and property damages.

**WARNING**

The replacement of lithium batteries or fuel cells by inadequately trained personnel could result in a hazard.

**WARNING**

- The battery has a safety device. Do not disassemble or modify the battery pack.
- Excessive vibration outside the final product, puncture, contact with metals, or tampering with the battery can cause it to fail.
- Do not short-circuit the battery pack by directly connecting the negative terminals with metal objects.
- Do not keep a battery pack in your pocket, purse, or anywhere together with other metal (conductive)objects.
- Never use force to install(insert) the battery pack.
- If the battery pack is damaged in any way stop using immediately.
- Immediately discontinue use of the battery if, while using, charging, or storing the battery, the battery emits an unusual smell, feels hot, changes color or shape, or appears abnormal in any other way.
- Only use GE HealthCare designated battery packs inside a GE HealthCare system and or GE HealthCare Charger, never use an aftermarket or non-GE HealthCare authorized battery pack, if you are unsure, contact GE HealthCare for confirmation.
- The battery contains safety and protection devices, which, if damaged, may cause the battery pack to fail.
- Do not heat the battery or discard it in a fire.
- Do not charge the battery near a heat source, such as a fire or heater.
- Do not leave the battery in direct sunlight.
- Do not pierce the battery with a sharp object, hit it, or step on it.
- Do not solder a battery.
- Do not connect the battery to an electrical power outlet.

Battery Storage

In order to ensure the capacity of the battery in the fully-charged state is not less than 80% of the initial fully-charged state. Need to store the battery at low temperature and low humidity, no dust and no corrosive gas atmosphere.

- Storage time $\leq$ 3 months: 0 - 45°C (32 - 113°F)
- Storage time $\leq$ 12 months: 0 - 23°C (32 - 73.4°F)
- Humidity: $<$ 65%



WARNING

- If the ultrasound system is not being used on a monthly basis, the battery needs to be removed during the lengthy non-use period.
- Store in a cool, dry place. Never store or charge batteries in extreme temperatures, such as a hot car; near heat sources or fire; in pressurized or very cold environments.



CAUTION

To avoid the battery bursting, igniting, or fumes from the battery causing equipment damage, observe the following precautions:

- Exposure to liquids can cause internal corrosion or damage to the cells or to the Battery Management System (BMS). The BMS protects the battery from overcharging, high self discharge, or imbalanced charging of the cells, any of which can present the possibility of hazards during charging or use.
- Do not immerse the battery in water or allow it to get wet.
- Do not put the battery into a microwave oven or pressurized container.
- If the battery leaks or emits an odor, remove it from all possible flammable sources.
- If the battery emits an odor or heat, is deformed or discolored, or in a way appears abnormal during use, recharging or storage, immediately remove it and stop using it. If you have any questions about the battery, consult GE HealthCare or your local representative.

Discharge/Charge Cycle

When the battery is stored for three months or more, the customer should perform one full discharge/charge cycle.

NOTE

A full discharge/charge cycle means the system is turned on using battery power until the battery loses its charge completely and the system shuts down. Plug the ultrasound system in until the battery is fully charged as indicated by a green LCD light.

Upon receipt of the ultrasound system and before first time usage, it is highly recommended that the customer perform one full discharge/charge cycle.

If the battery has not been used for >2 months, the customer is recommended to perform one full discharge/charge cycle. It is also recommended to store the battery in a shady and cool area with FCC (full current capacity).

Getting Started

One Full Discharge/Charge Cycle Process:

1. Full discharge of battery to let the ultrasound system automatically shut down.
2. Charge the ultrasound system to 100% FCC (full current capacity).
3. Discharge of ultrasound system for complete shut down (takes one hour for discharge).

When storing packs for more than 6 months, charge the pack at least once during the 6 month timeframe to prevent leakage and deterioration in performance.



WARNING

- After the battery is fully discharged, do not leave it fully discharged, charge battery as soon as possible.
- Only charge battery packs in authorized GE HealthCare devices or chargers, never use a third-party charger.
- Avoid unattended charging of battery packs, or where objects such as carpet, furniture, wood, vinyl floors, curtains or other flammable objects are nearby.

Safe Use

Follow below instruction to ensure safe use and adequate maintenance of rechargeable battery.

1. Regularly check the appearance of the cell to see if there is no abnormal deformation (bulge) etc.
2. Check the battery interface to see if there is no foreign body and deformation.
3. Regularly charge the battery if the system is not in use.


View current battery status

When the system is running on battery, there is a battery icon  displayed in the system status bar.

Total battery power remaining is shown on the battery icon. When there is no battery, the AC plug icon  is displayed in the system status bar.

Select the battery icon and the following information window appears:

Figure 1-3 Battery Status Message



a -Battery charging: - This field displays the current power source: AC power and Battery.

- When the battery is in use, this field displays the current power remaining capacity.
- When there is no battery, "**No Battery**" appears.
- When the battery is not in use, this field displays current capacity (**charging...**).

b -remaining %: - This field displays the current power remaining percentage.

When there is no battery, "**0**" appears.

If the battery is not in use, the battery icon  appears as being charged in the system status bar.

Battery power low warning

If the battery is in use and the battery power is lower than 40%, a warning message appears to warn the user that the battery power is low.

A red warning message "The battery is low, please charge your battery immediately" will display on the main screen.

When the estimated current power is lower than 30%, a warning will display on the screen to warn that the system is about to shut down.

NOTE

When the battery power is low and the user cannot charge the battery in time, the system automatically shuts down. This protects the whole system. You need to charge the battery by connecting to power source immediately before the system shuts down or you may lose useful information.

Gel Warmer (Option)

Gel Warmer is used for heating the gel to provide a comfortable scanning environment. There are 2 temperature level on the gel warmer.

Figure 1-4 Gel Warmer Temperature Option



When the setting is on Level I, the temperature is 28 degree; When the setting is on Level II, the temperature is 38 degree. The LED light will be on no matter the gel warmer is set on I or II.

NOTE

In order to protect patient safety, the highest temperature for gel should be less than 41 degree.

Wired Footswitch (Option)

You can attach the Footswitch to the system by connecting it to one of the USB port on the rear of the system.

**CAUTION**

To avoid damage of the cable, keep the cable away from the wheels. Disconnect the Footswitch before moving the system.

Figure 1-5 Footswitches



1. 1 Pedal Type Footswitch
2. Footswitch MKF 2-MED USB GP26

You can configure the function of Footswitch from the **Utility > Application > Settings > Footswitch**.

NOTE

For single footswitch, please use Middle to do setting.

**CAUTION**

When using the Footswitch, DO NOT hold down the footswitch pedal. Press and release the Footswitch pedal. Pushing and holding down the pedal behaves the same way as pushing and holding down a key on the keyboard.

Control Panel Map

Controls are grouped together by function for ease of use. See the callout for this figure.

Figure 1-6 Control panel map - example



- | | | |
|----------------------------|----------------------|-------------------------------|
| 1. Probe and Gel Holders | 10. PW Mode/Y | 18. Cursor |
| 2. Touch Panel | 11. PDI | 19. Clear key |
| 3. Power On/Off | 12. CF Mode/Z | 20. Measure key |
| 4. Menu knobs | 13. 3D/4D | 21. Body Pattern key |
| 5. User Configurable keys | 14. B Mode key | 22. Comment key |
| 6. Quad, Single, Dual keys | 15. Depth | 23. Print key |
| 7. Assistant key | 16. Ellipse/Zoom key | 24. Freeze key |
| 8. M Mode/X | 17. Whizz key | 25. Trackball, Trackball Keys |
| 9. CW key | | |

NOTE

The factory default settings can be modified in **Utility > System > User Configurable Key**.

Control panel adjustment



CAUTION

To avoid injury or damage, make sure nothing is within the range of motion before moving the control panel. This includes both objects and people.

The control panel position can be adjusted for easy viewing and ease of use.

To raise/lower the Control panel

1. Grab and hold the Up/Down lever under the control panel.

Figure 1-7 Up/Down lever



2. Pull up or push down forcefully to adjust the height of control panel.
3. Release the lever at the desired height.

NOTE

If the control panel is adjusted after a long period or for the first time, it may need more strength and produce some noise. It will become smoothly and easily after several times.

To swivel the Control panel

1. Grab the swivel lever under the control panel and hold it to swivel the control panel.
2. Release the lever at the desired position.

Indicator Light (For Versana Premier Lotus only)

Figure 1-8 Indicator Light



There is an indicator light on the control panel to indicate the system status.

- White -** System is in normal working status.
- Blue -** System is in standby status.
- Yellow -** System is in low battery status.

Physical A/N keyboard (Option)

Physical A/N keyboard is under the control panel. User can push the keyboard to project forward.

Figure 1-9 Push the Physical A/N keyboard



The standard alpha-numeric keyboard has some special functions.

Table 1-5 Special Key Function

Keyboard key	Function
Esc	Exit current display screen.
Help	Access Online help/user manual.

Keyboard key	Function
Arrow	Annotation Arrow.
Eject	Eject media.
Spooler	Activates DICOM Job Spooler screen.
Macro	Creates a Fast Key.
Macro	Plays a Fast Key.
Home	Move annotation cursor to home position; shift+key to set current annotation cursor position as the new home position.
Text Overlay	Switch between user text annotation overlays.
Grab Last	Activate the last selected data for edit.
Word Delete	Erase word associated with comment cursor.
Prt Scr	Print Screen.
Fn + Left/Right Keys	Adjust audio volume.
Fn + Up/Down Arrow Keys	Adjust brightness.

If you encounter a problem and cannot collect the logs immediately:

Table 1-6 Key for Collecting the Log

Keyboard key	Function
Alt+1 or Alt+2	Place a marker in the log.
Alt+D	Collect the logs.

Once the logs are collected, the service representative would be able to see the marker you added which will help to troubleshoot the problem.

Top/Sub Menu

The Top/Sub Menu contains exam function and mode/function specific controls.

Figure 1-10 Top/Sub Menu controls



The Top/Sub Menu contains adjustable knobs associated with it. The adjustable knobs are used to adjust the selected function. For example, to increase or decrease frequency, turn the knob clockwise or anticlockwise; the paddle switches are used to access and adjust the menu. The functionality of these controls change depending upon the currently displayed menu.

Touch Panel

The touch panel contains exam function, mode/function specific controls, digital keyboard, Digital TGC and Fifth Column for quick access.

Figure 1-11 Touch Panel




Table 1-7 Touch Panel controls

No	Name	Description
1	Patient	Press to enter new patient information.
2	Probe	Press to select appropriate probe.
3	Report	Press to generate report.
4	End Exam	Press to end current exam.
5	Digital TGC	Press to activate digital TGC.
6	Digital Keyboard	Press to activate digital keyboard.
7	Brightness	Press to adjust Display, Touch Panel, Control Panel brightness.
8	Sound Adjustment	Press to adjust volume.

No	Name	Description
9	User Defined Gesture	Press to view user defined gesture configuration. User can configure the gesture by going to Utility > System > User Defined Gesture. User can configure Screen Clean Mode under User Defined Gesture for cleaning the Touch Panel.
10	Utility	Press to enter Utility.
11	Fifth Column	Swipe the fifth column with one finger from right to left/left to right to change user configured preset page (Frequent Presets/Recent Presets/Customized Presets).

Digital Keyboard


User can use digital keyboard on the Touch Panel.

NOTE

Digital keyboard supports the same key combinations as the physical keyboard.

User can press the icon on the right upper corner of the Touch Panel to activate digital keyboard.


Figure 1-12 Activate digital keyboard



Getting Started

User can press the icon on the right bottom to exit digital keyboard.

Figure 1-13 Exit digital keyboard



Digital TGC

User can adjust TGC by activating digital TGC on touch panel, refer to *Figure 1-11 Touch Panel* on page 32.

Digital Longitudinal Gain Compensation (Digital TGC) amplifies the returned signal to correct for the attenuation produced when the beam penetrates tissue. Digital Longitudinal Gain Compensation (Digital TGC) balances the image so that the density of echoes is evenly distributed across the image. It enables the user to adjust B-mode gain compensation in the longitudinal direction.

Figure 1-14 Digital TGC




Table 1-8 Digital TGC controls

No	Name	Description
1	Standard Presets	Quick access to standard presets.
2	User Defined Presets	User can press and hold a button to save current user defined preset.
3	Touch TGC	Touch area for TGC adjustment.
4	Exit	Press to exit digital TGC.

Monitor

Monitor Display

Figure 1-15 Monitor Display Tour



Getting Started

1. Institution/Hospital Name, Date, Time, Operator Identification
2. Patient Name, Patient Identification
3. Power Output Readout
4. Probe Identifier. Exam Preset
5. Caps Lock: (lit when on), network connection indicator, system messages display, InSite status, InSite controls
6. Current date and time
7. Active Images screen
8. Delete Image
9. Save As Menu
10. Slide Show
11. Follow Up
12. Image Preview
13. Trackball Controls and Status
14. Image Clipboard
15. Measurement Window
16. Gray/Color Bar
17. Measurement Results Window
18. Measurement Calipers
19. Region of interest
20. Focal Zone Indicator
21. Imaging Parameters by Mode
22. Body Pattern
23. User Configurable Keys
24. P1, P2, P3 keys
25. Cine Gauge
26. Depth Scale
27. 6Tc-RS Probe temperature display
28. 6Tc-RS Probe angle display

To adjust the Monitor Position

The monitor can be swiveled and tilted.

Hold the bottom of the monitor when you adjust the position of the monitor.

Figure 1-16 Hold the bottom of the monitor



**CAUTION**

When flipping up the monitor from the flip down status, you can hold the upper corner of the monitor.


**CAUTION**

To avoid injury or damage, make sure nothing is within the range of motion before moving the monitor. This includes both objects and people.

Adjust Monitor by Software

1. Press **Cursor** key on the control panel. Move the trackball to position the cursor over the adjustment icon, then press **Set** to display the setting menu.

Figure 1-17 Adjust Monitor by Software



2. Adjust the Color Temperature, Contrast and Brightness by moving the trackball and using **Set** key.
3. Press **Exit** to exit the setting menu.

NOTE

Press **Reset** button can reset the three color temperature settings to default.

Moving the system

Before moving the system

When moving or transporting the system, follow the precautions below to ensure the maximum safety for personnel, the system, and other equipment.



CAUTION

When the system is not in use AND/OR before moving/transporting the system, make sure the monitor is flipped down to prevent system damage.



CAUTION

This product is not to be used on transport vehicles (e.g. ambulances, aircraft).



CAUTION

Handle carefully. A drop of more than 5 cm can cause mechanical damages.

1. Adjust the LCD monitor and control panel to their lowest positions. Flip down the LCD monitor before moving the system.



WARNING

Avoid injuring hands or fingers when adjusting the arm of monitor and adjusting the control panel.



WARNING

Avoid pinching hands and fingers during moving the system.

2. Press the Power On/Off switch to power off the system.
3. Unplug the power cord.
4. Wind the power cable around the cord hooks.

NOTE

To prevent damage to the Power Cord, DO NOT pull excessively on the cord or make sharp bends while wrapping.

5. All cables from off-board peripheral devices and the ethernet connection must be disconnected from the console.
6. Ensure that no loose items are left on the console.
7. Connect all probes to be used while off site. Ensure that probe cables are out of the way from the wheels and not protruding beyond the console. Use the probe management hooks located below the Control Panel to further secure the probe cables.
8. Store all other probes in their original cases or in soft cloth or foam to prevent damage.
9. Store sufficient gel and other essential accessories in the provided space.
10. Disconnect the footswitch from the console.
11. Unlock the wheels.

Wheels

Examine the wheels frequently for any obvious defects that could cause them to break or bind.

All the four wheels have independent brake pedals.

Figure 1-18 Wheel lock



1. Locked: Press the **ON** brake to lock the caster wheels.
2. Unlocked: Press the **OFF** brake to unlock the caster wheels.



MOVING HAZARD

Never move the system with locked wheels.



CAUTION

When two or more people are releasing the wheel controls with the front and rear wheels, take extra precaution to prevent unexpected movement which could result in possible toe injuries.



CAUTION

If you use/park the system on an incline, you **MUST** use the brakes on the wheel.

Moving the System

1. Always use the handles to move the system.

NOTE

User can choose front or rear handle to push the system for better sight and safety.

2. Take extra care when moving the system long distances and on inclines. Ask for help if necessary.

NOTE

Use two or more people to move the system on inclines or long distances. Failure to follow instructions could lead to possible injury and/or equipment damage.

Avoid ramps that are steeper than ten degrees to avoid tipping over the system. Utilize additional care and personnel when moving on a steep incline (>5 degrees) or loading the system into a vehicle for transport.

NOTE

Wheel chair ramps are usually less than five degrees.



CAUTION

DO NOT attempt to move the console using any cables or fixtures, such as the probe connectors.

3. Use the foot brake (pedal), located on the bottom of the system in the front, when necessary.
4. Do not let the system strike walls or door frames.
5. Use extra care when crossing door or elevator thresholds.
6. Once the destination is reached, lock the wheels.



CAUTION

The system weighs less than 45 kg (99.2 lbs.), without any probes or peripherals. To avoid possible injury and equipment damage:

- Be sure the pathway is clear.
- Limit movement to a slow careful walk.
- Use two or more persons to move the system on inclines or long distances.



CAUTION

DO NOT attempt to move the system with the monitor by pulling cables or belts placed around the monitor and/or monitor arm.



Do not put your whole body weight on the foot-rest holder.

Transporting the System

Use extra care when transporting the system using vehicles.

In addition to the instructions used when moving the system, also perform the following:

1. Before transporting, place the system in its special storage case.
2. Ensure that the system is firmly secured while inside the vehicle.
3. Secure the system with straps or as directed otherwise to prevent motion during transport.



CAUTION

Remove all peripherals from the console before transportation.

Load the system onto a vehicle

1. Only use vehicles that are designed for transport of the ultrasound system.
2. Load and unload the system to a vehicle parked on a level surface.
3. Ensure that the transporting vehicles can handle the weight of the system plus the passengers.
4. Ensure that the load capacity of the lift is capable of handling the weight of the system.
5. Ensure that the lift is in good working order.
6. Secure the system while it is on the lift so that it cannot roll. Use either wood chocks, restraining straps, or other similar types of constraints. Do not attempt to hold it in place by hand.
7. Employ two or three persons to load and unload safely from a vehicle.
8. Load the unit aboard the vehicle carefully and over its center of gravity. Keep the unit still and upright.

NOTE

Do not lay the unit down on its side.

9. Ensure that the system is firmly secured inside the vehicle.

NOTE

DO NOT restrain the Versana Premier/Versana Premier Lotus at the monitor or the monitor neck using a belt.

10. Prevent vibration damage by driving cautiously. Avoid unpaved roads, excessive speeds and erratic stops or starts.

System Start-Up

Connecting the System

To connect the system to the electrical supply:

1. Ensure that the wall outlet is the appropriate type.
2. Ensure that the power switch is turned off.
3. Unwrap the power cable. Make sure to allow sufficient slack in the cable so that the plug is not pulled out of the wall if the system is moved slightly.




CAUTION

Use the appropriate power cord provided by or designated by GE HealthCare.

4. Attach the power plug to the system, and then install the power cable bracket (b). Ensure the power plug is placed securely in power cable bracket.

Figure 1-19 Supply Mains Switch and Power Plug



5. Push the power plug securely into the wall outlet.

NOTE

Do not use an extension cord or adapter plug.

6. Switch ON the Supply Mains Switch (a). The On/Off button on the control panel illuminates blue, indicating that the system is in standby mode.



WARNING

To avoid risk of fire, the system power must be supplied from a separate, properly rated outlet.

Under no circumstances should the AC power plug be altered, changed, or adapted to a configuration rated less than specified. Never use an extension cord or adapter plug.

To help assure grounding reliability, connect to a "hospital grade" or "hospital only" grounded power outlet.



CAUTION

To ensure that the power cable does not disconnect during system use. If the system is accidentally unplugged, data may be lost.

**WARNING**

Failure to provide an adequate earth circuit can cause electrical shock, resulting in serious injury.

Connection of additional protective earth conductors or potential equalization conductors is not necessary in most cases and is only recommended for situations involving multiple equipment in a high-risk patient environment to provide assurance that all equipment is at the same potential and operates within acceptable leakage current limits. An example of a high-risk patient would be a special procedure where the patient has an accessible conductive path to the heart such as exposed cardiac pacing leads.


**CAUTION**

Only connect ultrasound system and mains-operated accessories to the appropriate wall outlet. DO NOT connect them to a single or multiple socket outlets, an extension cord, power strip or an adapter plug.

**WARNING**

To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

Figure 1-20 Example Plug and Outlet Configurations



1. 100-120 VAC, 10A Plug and Outlet Configuration
2. 220-240 VAC, 10A Plug and Outlet Configuration

NOTE

Country-specific power cords are currently available for Argentina, Australia/New Zealand, China, Denmark, India/South Africa, Switzerland, United Kingdom, Europe, USA, Israel, Brazil and Japan.

Voltage level check

Check the rating label on the rear side of the system. Check the voltage range indicated on the label.

Connecting to the electrical outlet

**WARNING**

POWER OUTAGE MAY OCCUR. The ultrasound unit requires a dedicated single branch circuit. To avoid circuit overload and possible loss of critical care equipment, make sure you DO NOT have other equipment operating on the same circuit.




CAUTION

Disconnect the power plug from the wall outlet or press the Power On/Off switch to turn the power off in case an emergency should occur. Ensure easy access to the power plug and the Power On/Off switch.

Supply Mains Switch

The Supply Mains Switch is located at the rear panel of the system. On supplies main power to all internal systems. Off removes main power from all internal systems.

Figure 1-21 Supply Mains Switch (located on the rear panel)



The Supply Mains Switch should stay in the **On** position; **DO NOT** hold the switch in the **On** position. If the Supply Mains Switch remains **On**, follow the Power **On** procedure.

1. On position
2. Off position

NOTE

If the Supply Mains Switch does not remain in the **On** position or trips again:

1. Disconnect the Power Cable.
2. Call Service immediately.

DO NOT attempt to use the system.

Power On



CAUTION

Press the Power On/Off switch to turn the power on. The Supply Mains Switch must also be in the on position. For Supply Mains Switch location, see *Supply Mains Switch* on page 44 for more information.

After the Supply Mains Switch is on, the LED light at the bottom of the monitor is lit and the Power On/Off switch illumination turns to white.

After a successful boot-up process, the Power On/Off switch illumination turns from white to green, the keyboard backlight is lit.

To turn on the system

1. Ensure that the unit is properly plugged into an AC outlet of sufficient capacity.
2. Turn on the breaker at the back of the system
3. Momentarily press the On/Off switch. After powered on successfully, the Power On/Off switch illumination is turned from white to green.
4. The system should now go through its boot-up process with no further user intervention.

Power Up Sequence

The system is initialized. During this time:

- The system boots up and the status is reflected on the monitor.
- Probes are initialized for immediate operation.
- Peripheral devices are activated on power up.



CAUTION

Please check that system date and time are correct after system is powered on.


After initialization is complete, controls on the Control Panel backlight and the default B-Mode screen is displayed on the monitor (if a probe is connected).

Login

Personal IDs and associated passwords can be preset on the ultrasound system.

If the User Auto Logon preset is blank, you are prompted to login.

Figure 1-22 Operator Login Window




1. **Operator:** Select the Operator.
2. **Password:** Enter Operator's password (optional).
3. Select **OK** or **Cancel**.
 - **OK:** Standard logon
 - **Cancel:** Cancel logon

NOTE

This ultrasound system **Local Disk** is encrypted by default. The user will be notified to complete encryption process by changing password and saving the recovery keys. If user did not complete this process, a warning message will appear each time the user logs in as an administrator (**ADM**).

Figure 1-23 Warning message for product encryption



Please navigate to **Utility > Admin > System Admin > Disk Encryption** to complete this process. For the detail information, please refer to **Basic Service Manual**.

4. **Reset Logon Accounts To Factory Default (Reboot needed):** Plug in SSA Dongle to reset Logon accounts.

If the user wants to skip the Login procedure, please check the **Use Auto Logon** preset from **Utility > Admin > Logon**.

Check System Date and Time

A warning message **Please check the system date and time are correct.** appears on the screen when the system is powered on. This warning message appears for the possible reasons:

- The system is not boot up for over 14 days.
- The system time has been changed by 24 hours earlier than the current system time of last boot-up.

A warning message pop out to remind the user to check the system date in case the system date and time is incorrect.

Move the cursor to **OK** and press Cursor key on the control panel to select **OK**. The system enters scanning mode.

Check the system date and time. If it is incorrect, follow below steps to reset the system date and time.

1. Enter **Utility > System > General > Date/Time**.
2. Reset the system date and time.

3. Select **Apply** and then select **OK**.
4. Select **Save**.

Power Off

To power off the system:

1. When you shutdown the system, enter the scan screen and lightly press the **Power On/Off** switch on control panel once. The System-Exit window is displayed.

NOTE

DO NOT press and hold down the Power On/Off switch to shutdown the system.

2. Using the Trackball, select **Shutdown**.

The shutdown process takes a few seconds and is completed when the Power On/Off switch illumination turns from green to white.

NOTE

DO NOT select **Exit** for shutdown. **Exit** is only available to service representative.

NOTE

If the system has not fully shut down in 60 seconds in the power-off sequence, press and hold down the On/Off switch until the system shuts down.

3. Disconnect the probes.

Clean or disinfect all probes as necessary. Store them in their shipping cases or another appropriate probe storage system to avoid damage.

NOTE

DO NOT turn off the Supply Mains Switch before the monitor display turns off.

Data may be lost or system software damaged if the Supply Mains Switch is turned off before the system shuts down.

4. Disconnect the power plug from the power outlet.

Log Off

To Logoff,

1. Press the On/Off switch momentarily and a **System Exit** window appears.
2. Press the Log Off switch to logoff current user account.

Press the On/Off Switch to fast boot up.

NOTE

If the user logoff the system, system will not save user configuration when log in next time.

Standby Mode

To activate Standby Mode:

1. Press the **Power On/Off** switch and select **Standby**.
2. A warning message display to ask user saving data before entering standby mode. Select **OK** to continue or select **Cancel** to return to system logon status.
3. The system will turn off all backlight and disable all keys (Function Keys, A/N keys, Trackball) except On/Off switch.

NOTE

The backlight of On/Off switch will flash slowly to indicate that system is in standby mode.

Press the On/Off Switch to fast boot up.



CAUTION

Low battery can cause system power off during Standby Mode. If the fast boot up failed, please plug in AC power and try again.

Crash Recovery Instructions

In cases where the system detects an internal error, the system may reboot on its own. If this happens, the system automatically returns to the start-up screens. All images and measurements, except for generic reports, are preserved in the system.

The system automatically ends the current exam and permanently store all the images and measurements. When the system reboots, check that all images and measurements have been preserved in the system. Then, simply hold down the power switch to initiate a normal power down sequence.

NOTE

If the image is not updated properly when the system is up, shut down the system again.

NOTE

Generic reports are not saved if the system crashes before you save it.



WARNING

A system crash may cause internal SSD (Solid State Disk) corruption. Avoid using the internal SSD as a permanent storage device. Backup data on a regular basis.

Database Protection Mechanism

Abnormal shutdown, such as power outage, has the risk to cause database corruption. Reboot the system after abnormal shutdown. The below dialog box displays when the system visits damaged database. Press **OK**, then the system restores the database which have been

backuper most recently in the normal condition. The data occurred since last time of normal store and normal shutdown may lost.

Probes

Connecting the Probe

Only use approved probes.



CAUTION

Inspect the probe before and after each use for damage or degradation to the housing, strain relief, lens, seal, cable and connector. DO NOT use a transducer which appears damaged until functional and safe performance is verified. A thorough inspection should be conducted during the cleaning process.



CAUTION

Remove any dust or foam rests from the probe pins.



CAUTION

Fault conditions can result in electric shock hazard. Do not touch the surface of probe connectors which are exposed when the probe is removed. DO NOT touch the patient when connecting or disconnecting a probe.

Probes can be connected at any time, regardless of whether the console is powered on or off.

To connect a probe:

1. Place the probe's carrying case on a stable surface and open the case.
2. Carefully remove the probe and unwrap the probe cord.
3. Put the probe in the probe holder.



CAUTION

DO NOT allow the probe head to hang free. Impact to the probe head could result in irreparable damage.

4. Hold the probe connector horizontally with the cable pointing upward.
5. Slide the connector lock to the left (unlocked position).
6. Insert the connector into the probe port and slide the connector lock to the right position to lock the probe connector.

NOTE

For the straight connector, please connect to the side port as illustrated on *Figure 1-24*.
Lock the probe connection on page 51.

Figure 1-24 Lock the probe connection



7. Wrap the probe cables around the probe cable hooks under the control panel. Ensure the probe cables are wrapped properly, not extended beyond sides of console.

Cable handling

Take the following precautions with probe cables:

- Keep free from wheels.
- Do not bend the cable acutely.
- Avoid crossing cables between probes.

Activating/Deactivating the Probe

To activate the probe, select the appropriate probe from the probe indicators in Probe screen.

To deactivate the probe, select another probe or press Freeze key to enter into freeze mode.

The probe's default settings for the mode and selected exam are used automatically.

**CAUTION**

Make sure that the probe and application names displayed on the screen correspond to the actual probe and application selection.

Disconnecting the Probe

Probes can be disconnected at any time. However, the probe should not be active when disconnecting the probe.

Getting Started

1. Deactivate the probe.
2. Slide the connector lock to the left position to unlock the probe.
3. Pull the probe connector straight out of the probe port carefully
4. Ensure the cable is free.
5. Be sure that the probe head is clean before placing the probe in its storage box.

Transporting Probes

Transport and store the ultrasound probes covered as needed and secured.

When transporting a clean or dirty ultrasound probe, ensure the probe is protected from cross contamination and possible damage. This can include the use of covers per the appropriate disinfection/contamination level, and the use of rigid containers or the probe holder on the ultrasound unit to secure.

NOTE

Do not use cloth or plastic bags to transport probes. This could result in damage to the probes.

When using a rigid transport case, such as the shipping container or a transport case, ensure that the probe is clean, and avoid damage to the probe by allowing nothing to protrude beyond the case when closing the lid, and secure the system connector in place so as not to damage the transducer head or lens.

- Secure the probe in its holder for moving short distances.
- When transporting a probe a long distance, store it in its carrying case.



WARNING

Placing an uncovered dirty or contaminated probe in a carrying case or shipping carton will contaminate the foam insert.



CAUTION

Avoid rapid and extreme temperature changes.

Storing the Probe

Prior to storage, it is essential to ensure the probe is completely clean and dry following disinfection.

If the probe is not immediately reused, store the probe in a manner that will protect and keep the probe from being re-contaminated. This may be accomplished by placing the probe lens upward on a wall mounted rack, probe holder on the ultrasound system, or in a storage cabinet with filtered air flow, using a disposable storage cover placed over the probe when needed. Avoid dangling the probe to prevent contact damage. For endocavitary probes, the optional horizontal probe holder mounted on the ultrasound machine can be used.

Do not use the shipping case or any closed container for other than short term transportation or shipping.



CAUTION

Avoid lengthy exposure of the probe to direct sunlight or to a strong ultraviolet light source.

It is recommended that all probes be stored in the provided carrying case for probe storage.

When storing a probe in the carrying case:

1. Place the probe connector into the carrying case.
2. Carefully wind the cable into the carrying case.
3. Carefully place the probe head into the carrying case. DO NOT use excessive force or impact the probe head.



CAUTION

DO NOT store probes in the side tray. To avoid damage, store the probe in its carrying case.

Preparing for an Exam

Patient Screen

User can go to **Utility > System > General > Patient Workflow (restarted needed)** to select Traditional or Simplified patient workflow.

Press the **Patient** key on the control panel to display the Patient Screen on the monitor.

- Enter Patient Data with the alphanumeric keyboard.
- To navigate through the Patient Entry menu, use the **Tab** key or **Trackball** and **Set** to move and fix the cursor.

Patient Screen - Traditional

Figure 1-25 Patient Screen - Traditional

Table 1-9 Patient Screen - Traditional

No.	Function	Description
1.	Image Management	• Patient - Provides a search and creation of patient. (currently selected). • Image History - Provides a list of images per exam for the currently selected patient. • Active Images - Provides image preview of the currently selected exam. • Data Transfer - Provides an interface to handle patient data from a remote device.

No.	Function	Description
2.	Function Selection	- New Patient - Used to clear patient entry screen in order to input a new patient's data into the database. - Register - Used to enter new patient information into the database prior to the exam. NOTE If you are using the auto-generate Patient ID feature, do not select Register. It is always a good practice to Register all patients. - Details - Select the Details box to activate/deactivate the exam details. The Exam Description pull-down selection is used as the DICOM identifier.
3.	EZBackup/EZMove	- EZBackup - One-step method to backup patient images to an external media. - EZMove - Move and delete patient images.
4.	Dataflow Selection	Select archive and other pre-defined services. If you use a DVD-R, select DICOM CD Read in Dataflow. If you place the cursor on the icon, the pop-up menu displays disk capacity.
5.	Scan	Used to exit Patient Menu and enter scanning screen.
6.	Patient Information	- Patient ID Number (User cannot re-edit ID after once edition.) - Other ID Number. The system now allows you to enter a second identification number for the same patient, which may be required in certain countries. This is only displayed if enabled on the Utility > Connectivity > Miscellaneous screen. - Patient Name-Last, First and Middle - DOB (Birthdate) - Age (automatically calculated when birthdate is input) - Gender
7.	Category Selection	Select from 10 exam application categories. When a category is selected, the measurement and category presets are displayed.
8.	Exam Information	Shows the Current/Active Exam information. Information pertinent to the selected exam category appears in the window. All possible information needs to be entered. - Clear-Clears existing data. - Past Exam (only for OB)-Input past exam data (register the patient before using).
9.	Patient View/Exam View	Display either patient or examination list.

Patient View

Figure 1-26 Patient View




Table 1-10 Patient View

No.	Function	Description
1.	Patient View	List the patients in the database. When you double-click the patient on the patient list using the Set key, the Review screen or New Exam entry screen displays depending on the preset. You set this preset via Select Review or New Exam in Utility > Connectivity > Miscellaneous > Double click on patient list to start preset.
2.	Search key	Select a search criteria. NOTE Criteria "Img. Archived" means that the exam was backed up to external media by EZBackup or Export.
3.	String	Enter appropriate information for search criteria. NOTE If you select Locked (Y, N) or Archived (Y, N) for the Search key, enter Y (Yes) or N (No).
4.	Clear	Clears the entered string.
5.	Hide	Hide the patient view.
6.	Review	Highlight the patient and press Review. Exam view is displayed.. NOTE If the selected patient has a Current Exam or the selected exam is the Current Exam, the Review button changes to "Resume Exam" on Patient List.

No.	Function	Description
7.	New Exam	Creates a new exam based on a current or searched patient.
8.	Delete	Deletes Patient/Exam. NOTE "Delete" is only displayed when you login as Administrator.
9.	Lock/Unlock	Locks the exam/patient. Prevents move and delete functions. If you select the patient, all exams are locked. If you select one exam, the selected exam is locked and the lock icon displays in the patient ID cell.
10.	Listing XX of XXX	Displays the quantity of patients in the search window and the quantity of patients in the database.

Exam View

Figure 1-27 Exam View




Table 1-11 Exam View

No.	Function	Description
1.	Exam View	The system can display the Detail Mode instead of Exam View when you double click the patient in the patient list or select the patient in the patient list and then select Details. If the Detail Mode preset on Utility > Connectivity > Miscellaneous menu is selected, the Detail Mode displays.
2.	Disk	Displays the disk name on which you saved the exam's image data. If "+" displays behind the disk name, the data is saved on two or more disks.
3.	Delete Exam	Delete one or more exams from Exam View. NOTE "Delete" is only displayed when you login as Administrator.
4.	Send To	Send the images to the DICOM Device.

Patient Screen - Simplified

Figure 1-28 Patient Screen - Simplified

The screenshot displays the 'Patient Screen - Simplified' interface. It features a top navigation bar with 'Patient', 'Archive', and 'Image Review' options. The main area is divided into several sections:

- 1**: A sidebar menu on the left with icons for Patient, Exam, and Report.
- 2**: A 'Scan' button at the bottom left.
- 3**: The 'Patient Info' section containing fields for Patient ID (230821-105800), Last Name, First Name, Middle Name, Gender (Male/Female), Age, and DOB.
- 4**: The 'Exam' section with tabs for ABD, OB, GYN, CAUD, VAS, UR, SMP, PED, MUSK, and THOR.
- 5**: A dropdown menu for 'Exam Description'.
- 6**: The 'Probe/Preset' section with buttons for various presets like Abd, Abd_Pen, Renal, Aorta, OB23, OB1, Fetal, and others.
- 7**: An 'Exam List' table at the bottom with columns for Date, Category, Exam Description, Img. size, M&A, and Disk.

Table 1-12 Patient Screen - Simplified

No.	Function	Description
1	Patient	- New Patient - Provides creation of patient. - New Exam - Creates a new exam based on a current or searched patient. - Detail - Used to check the detail of patient information and exam information. - Report - Used to enter report of current exam.
2	Function Selection	- Scan - Used to exit Patient Menu and enter scanning screen. - End Exam - Used to end current patient.

No.	Function	Description
3	Patient Information	- Patient ID: (User cannot re-edit ID after once edition.) - Default ID(Auto ID or Manual ID): It is generated by system automatically or inputted manually and only can be re-edited in Patient Info page once. - Temp ID(Simplified Patient only): It can be re-edited once only after saving images without patient ID etc. information, i.e. in emergency. The ID is shown 'xxxTEMP' that differs from Default ID. - Other ID: The system now allows you to enter a second identification number for the same patient, which may be required in certain countries. This is only displayed if enabled on the Utility > Connectivity > Miscellaneous screen. - Patient Name-Last, First and Middle - DOB (Birthdate) - Age (automatically calculated when birthdate is input) - Gender
4	Category Selection	Select from 9 exam application categories. When a category is selected, the measurement and category presets are displayed.
5	Exam Information	Shows the Current/Active Exam information. Information pertinent to the selected exam category appears in the window. All possible information needs to be entered. - Clear-Clears existing data. - Past Exam (only for OB)-Input past exam data (register the patient before using).
6	Probe/Preset	Display or select probe and preset information.
7	Exam List	Display examination list.

Archive

Figure 1-29 Archive




Table 1-13 Archive

No.	Function	Description
1.	Archive	List the patients in the database. When you double-click the patient on the patient list using the Set key, the Review screen or New Exam entry screen displays depending on the preset. You set this preset via Select Review or New Exam in Utility > Connectivity > Miscellaneous > Double click on patient list to start preset.
2.	Search key	Select a search criteria. NOTE Criteria "Img. Archived" means that the exam was backed up to external media by EZBackup or Export.
3.	String	Enter appropriate information for search criteria. NOTE If you select Locked (Y, N) or Archived (Y, N) for the Search key, enter Y (Yes) or N (No).
4.	Clear	Clears the entered string.
5.	Dataflow Selection	Select archive and other pre-defined services. If you use a DVD-R, select DICOM CD Read in Dataflow. If you place the cursor on the icon, the pop-up menu displays disk capacity.
6.	Add Exam	Creates a new exam based on a current or searched patient.
7.	Resume	Continues the exam for that patient if you select the last exam of the day.

No.	Function	Description
8.	EZBackup/EZMove	- EZBackup - One-step method to backup patient images to an external media. - EZMove - Move and delete patient images.
9.	Printer	Print the search list to a standard printer. Highlight the patient and press Cursor key. Select Print from the pop-up and press the right Set key.
10.	Send To	Send one or multiple patients' images to the peripherals in JPG/AVI, JPG/WMV, DICOM, Raw Dicom or Mpeg Vue format. DICOM - Digital Imaging and Communications in Medicine is a standard protocol for the management and transmission of medical images and related data. Raw DICOM - Based on DICOM images and related data, Raw DICOM contains GE proprietary information which can be post-processed on GE system.
11.	Import/Export	Select to import or export patient information (including Report, images and other patient information).
12.	Listing XX of XXX	Displays the quantity of patients in the search window and the quantity of patients in the database.
13.	Q/R	Allows for One Click Q/R (Query/Retrieve) for the current patient.
14.	Delete	Delete one or more exams from Exam View. NOTE "Delete" is only displayed when you login as Administrator.

Image Review

Figure 1-30 Image Review




Table 1-14 Image Review

No.	Function	Description
1.	Image Review	Select to review all images of the selected patient exam.
2.	Active Image	Provides preview of the currently selected exam.
3.	Image History	Provides a list of images per exam for the currently selected patient.
4.	Review	Highlight the image and press Review.
5.	Compare	Enter Follow-up function.
6.	Printer	Print the search list to a standard printer. Highlight the patient and press Cursor key. Select Print from the pop-up and press the right Set key.
7.	Send To	Send one or multiple patients' images to the peripherals in JPG/AVI, JPG/WMV, DICOM, Raw DICOM or Mpeg Vue format. DICOM - Digital Imaging and Communications in Medicine is a standard protocol for the management and transmission of medical images and related data. Raw DICOM - Based on DICOM images and related data, Raw DICOM contains GE proprietary information which can be post-processed on GE system.
8.	Delete	Delete one or more exams from Exam View. NOTE "Delete" is only displayed when you login as Administrator.

Scanning a New Patient

Warnings and Cautions



WARNING

Imaging functions may be lost without warning. Develop emergency procedures to prepare for such an occurrence. Failure to prepare for unexpected loss of functionality could lead to patient injury.



WARNING

Always ensure you have selected a dataflow. If No Archive is selected, no patient data is saved and a rescan may be required. A Ø appears next to Dataflow if No Archive is selected.



CAUTION

To avoid patient identification errors, always verify the identification with the patient. Make sure the correct patient identification appears on all screens and hard copy prints.

**CAUTION**

Always use the minimum power required to obtain acceptable images in accordance with applicable guidelines and policies.

**CAUTION**

Always use the system on a flat surface in the patient environment.

**CAUTION**

Ensure that the hands of the patient are away from the system during the exam.
The position of the operator and the patient vary by the anatomy of interest being scanned.
In most cases, the operator sits/stands straight in front of the operator console and the patient lies on the bed on the right (or left) side of the system.

**CAUTION**

Before each use of intracavitary probe, or after changing its observation mode/setting, please check to ensure that the image is a live (not the stored image) and has the correct image orientation.

**CAUTION**

The intracavitary probes supported by this product are only expected to contact intact mucosa.

Creating a new patient record

When starting a new patient's exam, ensure you do the following:

1. Select **Patient** on the touch panel.
2. Select **New Patient** on the patient menu.
If there are images on the clipboard, a pop-up menu appears. Specify whether you want to store images, delete images, or go to active images.
3. Choose the exam category.
4. Verify the dataflow.

NOTE

DO NOT use the removable media Dataflows on the New patient menu.

NOTE

The system will display a warning dialog when the patient is registered to "No Archive". If the "Warn register to No Archive" preset is selected in the **Utility > Connectivity > Miscellaneous** menu, a warning displays. Please select a different dataflow for permanent storage of patient data.

5. Fill in patient information.

NOTE

You can also select a patient from the patient database at the bottom of the Patient menu if the patient has a patient ID.

NOTE

Columns drive the ordering of the patients displayed. The columns that you select drive the order of the displayed patient database.

NOTE

Do not use the following characters when filling in patient information: " ` / : ; . , * < > | + = [] & ! , @ , # , \$, % , ^ , & , * , (,) , ? , / , ~ , [,] , { , } .

6. Select **Register**. Enter Past OB Exam information, if desired.
7. Select the probe to start scanning (or select Esc, Scan, or Freeze).
8. Perform the exam.

Perform an exam

1. Select the probe, exam category and application.
2. Perform an exam.
3. Store the raw data to the clipboard.

To store the still image, press **Freeze** and run the cineloop using the **Trackball**. Select the frame and press **Store**.

To store the cineloop, press **Freeze** and run the cineloop using the **Trackball**. Select the start/end frame and run the selected loop. Press **Store**.

Ending an exam

1. When you have completed the study, press **End Exam**.
2. The image management screen displays. Select the images (still frame or cineloop) you want to store or select **Select All** to store all images. Select **Permanent Store** to store the images permanently.



CAUTION

After completing the measurement, verify that the measurement result window is updated before you send or save the image.



CAUTION

DO NOT use the following special characters when saving images: ! , @ , # , \$, % , ^ , & , * , (,) , | , ; , : , < , > , ? , / , ~ , [,] , { , } .

Entering a Patient List

All patient information can be entered before starting an exam.

1. Press **Patient** to display the Patient Screen.
2. Press **New Patient** to erase the current patient data.
3. Select **Store All** if data from previous patient was not saved.
4. Enter the Patient ID.
5. Enter the patient and exam information.

6. Press **Register**.
7. Repeat above steps as required.

Select the patient from the Patient List. Select **Resume Exam** to continue the last exam that was performed on the selected patient.

Select **New Exam** to start a new exam on the selected patient.

Starting a new exam on an existing patient

1. Select **Patient** on the control panel.
2. Select the patient from the Patient List.
3. Select **New Exam**.
4. A new exam is created. Enter the data and begin the scan.

Scanning without entering any patient data

To scan a patient without entering any patient data until the end of the exam:

1. Scan the patient and save images to clipboard without the patient information, the system displays a Warning message, "A patient must be selected for permanent storage of images." Select **OK**, a warning message is also displayed at the bottom of the image monitor in red.
2. When the scanning is finished, press **Patient** to display the Patient Search screen.
3. Enter the Patient ID, patient data and exam information as necessary, select **Register**.
4. If images or measurements have been stored to the clipboard, the system will display the following message:
"Unsaved images, measurements or fetus number will be linked to the current patient information, continue?" Press **OK** if you want to permanently store the images/measurements that were just taken.
5. Enter the Active Images Screen, select **Permanent Store**.
6. Return to patient page, select **New Patient**.

When the system shuts down during an exam, the unsaved images and data are stored in temporary buffer. Next time the system starts up, the temporary buffered data is still available.

Ending a Patient

To end a patient:


Press **Patient > New Patient**, then select **Store All** in the pop-up menus to store exam data.

Changing Current Patient to Existing Patient (with Patient ID)

To change the current patient (with patient ID) to an existing patient, when there are some unsaved images on the clipboard for the current patient:

Select the existing patient from the patient list, the following dialog displays.

Figure 1-31 Unsaved Exam Data




1. **Store All** - All the unsaved images will be saved to the current patient.
2. **Review Images** - Review the unsaved images in Active Images page and select to permanent store or delete.
3. **Delete All** - Delete all the unsaved images.
The system will display a confirmation message: "Do you really want to delete all temporary images?" Select **OK** to delete the unsaved images, select **Cancel** not to delete the images.
4. **Cancel** - Do nothing with the unsaved data, and return to Patient page.

Changing Current Patient to Existing Patient (without Patient ID)

To change the current patient (with patient ID) to an existing patient, when there are some unsaved images on the clipboard for the current patient:

Select the existing patient from the patient list, the following dialog displays.

Figure 1-32 Unsaved Exam Data



1. **Store All** - All the unsaved images will be saved to the current patient. The current patient does not have Patient ID, the system will indicate to enter patient ID.
 - **Input manually** - Input the Patient ID and other information manually, and then select Register to register the new patient. And the system will display:
"Unsaved images, measurements or fetus number will be linked to the current patient information, continue?" Select OK to save the images.
 - **Auto Generate** - The system will generate a new patient ID and the unsaved images will be saved to this new patient automatically.
2. **Review Images** - Review the unsaved images in **Active Images** page and select to permanent store or delete.
3. **Delete All** - Delete all the unsaved images.
The system will display a confirmation message: "Do you really want to delete all temporary images?" Select **OK** to delete the unsaved images, select **Cancel** not to delete the images.
4. **Cancel** - Do nothing with the unsaved data, and return to **Patient** page.

Delete the existing patient/exam/image



CAUTION

Before deleting a patient or image from the Patient Screen, make sure you have already saved the data with EZBackup/EZMove, Backup, or Export. Verify the media before deletion.

Delete the existing patient

1. Search and select the patient in the patient list.
2. Select **Delete**. The confirmation dialog box displays.
OR
Press the **Cursor** key. A pop-up menu displays. Select **Delete**. The confirmation dialog box displays.
3. Select **OK** to delete or **Cancel**.

Delete multiple patients from the patient list

1. Select the multiple patients to be deleted from the patient list by holding down the Control key and selecting patients one by one or holding down the Shift key to select a group of patients.
2. Select **Delete**. The confirmation dialog box displays.
OR
Press the **Cursor** key. A pop-up menu displays. Select **Delete**. The confirmation dialog box displays.
3. Select **OK** to delete or **Cancel**.

Delete the existing exam

1. Search and select the patient in the patient list under patient screen.
2. Select **Exam View** to display patient exam screen.
3. Select the exam to be deleted.
4. Select **Delete Exam**. The confirmation dialog box displays.
5. Select **OK** to delete or **Cancel**.

Delete the existing image

1. Search and select the patient in the patient list under patient screen.
2. Select **Exam View**. The patient exam screen displays.
3. Select the exam which contains the image to be deleted.
4. Select **Active Images** to display the image list.
5. Select the image to delete and select **Delete**. The confirmation dialog box displays.
6. Select **Yes** to delete or **No** to cancel.

Barcode Reader

The Barcode Reader supports 1D and 2D encoding symbologies. The protocols supported are shown as below:

Table 1-15 Protocols supported by Barcode Reader

Category	Protocol
1D	Code 128, Code 93, Interleaved 2 of 5
2D	PDF417, Data Matrix

Activate Barcode Reader

The barcode reader needs to be activated before the first use.

To activate the barcode reader, follow the steps below:

1. Power on the system and connect the barcode reader to USB port on the ultrasound system.

NOTE

The following barcodes **MUST** be printed on paper before being scanned by the barcode reader.

2. Scan the four barcodes.
First, scan the Remove Custom Defaults barcode (1).
Second, scan the Activate Defaults barcode (2).
Third, scan the To enable 1900 PID code(3).

Finally, scan the USB Serial barcode (4).

Figure 1-33 Barcode for scanning



Setup Barcode Reader

Before using barcode reader, follow below instruction to complete configuration.

Select **Utility** > **Barcode** to enter configuration page.

Input Mode

There are **Patient ID**, **Complexation** and **Off** mode listed under Input Mode.

Off

If **Off** is selected, user needs to use the keyboard to enter the Patient ID.

Patient ID

If **Off** is selected, scan the Barcode as the Patient ID or enter the Patient ID using the keyboard..

Complexation

If Complexation is selected, user can scan the barcode to input the Patient demographics or use the keyboard to enter the Patient demographics.

To enter the Patient demographics:

1. Enter a barcode string in the **Input Data** tab by scanning from a barcode or typing with the keyboard.

The following items can be included in the barcode:

- Patient ID
- First Name, Last Name, Middle Name
- Birth Year, Birth Month, Birth Day

NOTE

The length of Birth Year contains 4 digits, Birth Month contains 2 digits and Birth Day contains 2 characters; they should always be provided together.

- Gender

NOTE

Please input **F** for female and use **M** for male.

2. Configure the Start and End position for each item.

NOTE

If the barcode does not contain the information of any item, configure the Start and End position as "0".

For example, if the scanned barcode is "0102LastNameFirstNameMiddleName20001212F", the configuration and results display as the following:

Figure 1-34 Barcode Manager

The screenshot shows the Barcode Manager interface. At the top, there are radio buttons for 'Patient ID', 'Complexation', and 'Off'. Below that, there are radio buttons for 'Patient ID', 'Accession #', and 'Off'. A purple box highlights the barcode string '0102LastNameFirstNameMiddleName20001212F'. Below the barcode string, there is a dropdown menu for 'Barcode Scanner' set to 'Xenon 1900 Area-Imaging Scanner'. The 'Input Data' field contains the same barcode string. Below this, there are several input fields for patient information, each with a dropdown menu for the start and end positions. The fields are: Patient ID (1, 2, 01), Birth Year (32, 35, 2000), Other ID (3, 4, 02), Birth Month (36, 37, 12), Last Name (5, 12, LastName), Birth Day (38, 39, 12), First Name (13, 21, FirstName), Gender (40, 40, F), Middle Name (22, 31, MiddleName), Male (M), and Female (F). At the bottom, there are 'Save', 'Exit', and 'Cancel' buttons.

Table 1-16 Complexation mode Input Data configuration

Item	Start Position	End Position	Input Data Element
Patient ID	1	2	01
Other ID	3	4	02
Last Name	5	12	Last Name
First Name	13	21	First Name
Middle Name	22	31	Middle Name
Birth Year	32	35	2000
Birth Month	36	37	12
Birth Day	38	39	12
Gender	40	40	F

- Click **Save** to save the configuration.

Chapter 2

Performing an Exam

In this section

Optimizing the Image.	74
Measurement and Analysis.	95

Optimizing the Image

Optimizing B-Mode

B-Mode is intended to provide two-dimensional images and measurement capabilities concerning the capabilities of anatomical structures.

Figure 2-1 B-Mode Display -- Representative Example



Typical B-Mode Exam Protocol

A typical examination using B-Mode might proceed as follows:

1. Record exam-related patient information. Verify system setup (probes and presets).
2. Position the patient and the console for optimum operator and patient comfort. Perform the scan.
3. Complete the study by collecting all the data.

B Mode Controls

Table 2-1 B Mode Controls

Control	Bioeffect	Description/Benefit
Depth	Yes	Depth controls the distance over which the B-Mode images anatomy, and the field of view. To visualize deeper structures, increase the depth. To visualize the flatter structures, decrease the depth.
Gain	No	B-Mode Gain increases or decreases the amount of echo information displayed in an image. It may have the effect of brightening or darkening the image if sufficient echo information is generated.

Control	Bioeffect	Description/Benefit
Focus	Yes	Increases the number of focal zones, moves the focal zone(s) and change the zone width so that you can tighten up the beam for a specific area. A graphic caret corresponding to the focal zone position(s) appears on the right edge of the image.
Whizz	No	Whizz in B mode (ATO: Auto Tissue Optimization) will continuously optimize the brightness, contrast and uniformity of B mode images when scanning on different tissues. Whizz in PW/CW Doppler Mode (ASO: Auto Spectral Optimization) optimizes the spectral data. Auto adjusts the Velocity Scale/PRF (live imaging only), baseline shift, and invert (if preset). Upon deactivation, the spectrum is still optimized. Specify the Whizz Level: Low or High via Utility > Whizz/ Application > Whizz > Whizz Level.
CrossXBeam	Yes	CrossXBeam is the process of combining three or more frames from different steering angles into a single frame. CrossXBeam is available on Convex and Linear probes. CrossXBeam combines multiple co-planar images from different view angles into a single image at real-time frame rates, using bi-cubic interpolation.
Coded Harmonic Imaging (Phase Inversion Harmonic) (CHI)	Yes	Harmonic imaging utilizes Digitally Encoded Ultrasound (DEU). Coded Harmonics enhances near field resolution for improved small parts imaging as well as far field penetration.
Frequency	Yes	Adjust Frequency until the desired frequency is selected. Multi Frequency mode lets you downshift to the probe's next lower frequency or shift up to a higher frequency.
Steer	Yes	You can slant the B-Mode or Color Flow acoustic beam without moving the probe. The steer function only applies to linear probes. NOTE: Steer is not available while using CrossXBeam.
Virtual Convex	Yes	Virtual Convex provides a larger field of view in the far field. The Virtual Convex feature provides an extra wide field of view on Linear and sector probes, by introducing an offset, which is similar to radius of convex, and it defines the FOV of Linear and Sector probes. Virtual Convex allows for a wider field of view. Available in B-Mode, Color Flow Mode, and Doppler Mode. CrossXBeam is available on Virtual Convex with linear probes. Activating Virtual Convex may change the TI and/or MI. Observe the output display for possible effects.

Control	Bioeffect	Description/Benefit
TGC	No	TGC amplifies returning signals to correct for the attenuation caused by tissues at increasing depths. TGC slide pots are spaced proportionately to the depth. The area each pot amplifies varies as well. A TGC curve may appear on the display (if preset), matching the controls that you have set (except during zoom). You can choose to deactivate the TGC curve on the image.
Width	Yes	You can widen or narrow the size of the sector/convex angle to maximize the image's region of interest (ROI).
Tilt	Yes	You can steer the sector angle to get more information without moving the probe while in B-Mode, M-Mode, Doppler Mode, and Color Flow Mode. NOTE Tilt is not available on linear probe. NOTE Tilt is available when CrossXBeam is off.
Revert	No	Flips the image 180 degrees left/right.
Dynamic Range	No	Dynamic Range controls how echo intensities are converted to shades of gray, thereby increasing the adjustable range of contrast.
Line Density	Yes	Optimizes B-Mode frame rate or spatial resolution for the best possible image.
Gray Map	No	The system supplies B, M, and Doppler Mode system maps.
Frame Average	No	Temporal filter that averages frames together, thereby using more pixels to make up one image. This has the effect of presenting a smoother, softer image.
Colorize	No	Colorize is the colorization of a conventional B-Mode image or Doppler Spectrum to enhance the user's ability to discern B, M, and Doppler Mode intensity valuations. Colorize is NOT a Doppler Mode. NOTE You can colorize real-time or CINE images or Timeline CINE, but not DVR images. Colorizes the gray scale image to enhance the eye's discrimination capability. The colorize bar displays while Colorize is activated.
Edge Enhance	No	Edge Enhance brings out subtle tissue differences and boundaries by enhancing the gray scale differences corresponding to the edges of structures.

Control	Bioeffect	Description/Benefit
LGC	No	LGC (Lateral Gain Compensation) balances the image so that the density of echoes is the same throughout the image. It enables user to adjust B mode gain compensation along lateral direction.
Speed of Sound	No	Speed of Sound is available on all probes for all applications. The Speed of Sound control can be configured via Utility > Application > Imaging Controls. With SoS adjustment, user can adjust the SOS close to the real ultrasound transmission speed in different tissues and patients, which provides better focusing to optimize the image quality. NOTE Speed of Sound displays on the monitor display as "SoS" with the speed following, "SoS 1500" (when the speed of sound is not equal to default value 1540 m/s).
Rotation	No	You can flip the image 90 degrees/180 degrees/270 degrees. CAUTION: When reading an rotated image, be careful to observe the probe orientation to avoid possible confusion over scan direction or left/right image reversal.
Rejection	No	Selects a level below which echoes will not be amplified (an echo must have a certain minimum amplitude before it will be processed).
Suppression	No	Suppresses the noise in the image.
SRI-HD	No	SRI-HD (High Detection Speckle Reduction Imaging) is an adaptive algorithm to reduce the unwanted effects of speckle in the ultrasound image. Image speckle usually appears as a grainy texture in otherwise uniform areas of tissue. Its appearance is related to image system characteristics, rather than tissue characteristics, so that changes in system settings, such as probe type, frequency, depth, and others, can change the appearance of the speckle. Too much speckle can impair image quality and make it difficult to see the desired detail in the image. Likewise, too much filtering of speckle can mask or obscure desired image detail. Extra care must be taken to select the optimal SRI-HD level. Tips and Notes When selecting the SRI-HD level, observe the effect of SRI-HD in the desired region of interest and make a real-time comparison with the original image. The optimal level depends on the clinical situation. Observing the original and SRI-HD-processed images together helps to determine whether too much or too little SRI-HD has been applied. Dual image mode for SRI-HD can also be activated on a stored CINE Loop. This allows you to always see the original, unprocessed or enhanced image by going into the Dual display mode and to change the SRI-HD settings when reviewing the CINE Loop.

Control	Bioeffect	Description/Benefit
LOGIQView (Option)	No	LOGIQView provides the ability to construct and view a static 2D image which is wider than the field of view of a given transducer. This feature allows viewing and measurements of anatomy that is larger than what would fit in a single image. Examples include scanning of vascular structures and connective tissues in the arms and legs. LOGIQView constructs the extended image from individual image frames as the operator slides the transducer along the surface of the skin in the direction of the scan plane. The quality of the resulting image is somewhat user-dependent and requires some additional skill and practice to develop proper technique and become fully proficient. LOGIQView is only available in B mode.
Full Screen	No	Full Screen button is used to enlarge image and black area on the monitor. Only Full Screen button and B button can be used in Full Screen mode.

Digital LGC (Lateral Gain Compensation)

Digital LGC (Lateral Gain Compensation) amplifies returning signals to correct for the attenuation caused by tissues when the beam penetrates.

Digital LGC (Lateral Gain Compensation) balances the image so that the density of echoes is the same throughout the image. It enables user to adjust B mode gain compensation along lateral direction.


Adjusting LGC

User can activate LGC from the Touch Panel.

Figure 2-2 LGC icon on the Touch Panel




Figure 2-3 LGC ON



To adjust and save User Defined LGC:

Figure 2-4 Adjusting the LGC



See below description for the usage of the icons marked on Figure 2-4:

- a.** Adjust the slide pots to your preferred value.
- b.** Press to reset the LGC value to the center position.
- c.** User Defined LGC settings.

NOTE

Press and hold a button to save user defined preset.

- d. Press to exit LGC.

Optimizing M-Mode

M-Mode is intended to provide a display format and measurement capability that represents tissue displacement (motion) occurring over time along a single vector.

Typical exam protocol

A typical examination using M-Mode might proceed as follows:

1. Get a good B-Mode image.
2. Press **M-Mode**. Survey the anatomy and place the area of interest near the center of the B-Mode image.
3. Move the trackball to position the mode cursor over the area that you want to display in M-Mode.
4. Adjust the Sweep Speed, TGC, Gain, Power Output, and Focus Position, as needed.
5. Press **Freeze** to stop the M trace.
6. Record the trace to disk or to the hard copy device.
7. Press **Freeze** to continue imaging.
8. To exit, press M-Mode.

M-Mode Controls

NOTE

You can set default value of each parameter by probe and application on the **Utility > Imaging** page.

Table 2-2 M-Mode Controls

Control	Bioeffect	Description/Benefit
Sweep Speed	No	Changes the speed at which the timeline is swept. Available in M-Mode, Doppler Mode and M Color Flow Mode.
Anatomical M-Mode (option)	Yes	Anatomical M-Mode gives you the ability to manipulate the cursor at different angles and positions. The M-Mode display changes according to a motion of the M cursor.

Anatomical M-Mode (AMM) and Anatomical Color M-Mode (ACMM)

Anatomical M-Mode allows manipulation of the cursor at different angles and positions. The M-Mode display changes according to the position of the cursor.

Anatomical M-Mode displays a distance/time plot from a cursor line, which is independent from the axial plane. AMM is available in B, Color and TVI.


NOTE

To set up AMM, go to Utility > Imaging > AMM. Select the specific probe and parameters.

NOTE

Turn off CrossXBeam before activating AMM/CAMM.

Figure 2-5 Anatomical M Mode

**Activating**

To activate Anatomical M-Mode, while in M-Mode, press the **Anatomical M** Touch Panel control.

NOTE

AMM is not available on Linear probes.

NOTE

Anatomical M-Mode can also be used with previously acquired digitally stored 2D images. More than one heart cycle should be stored if performing M-Mode in post processing.

To activate Anatomical Color M-Mode, after accessing Anatomical M-Mode, activate Color Flow (CF).

Adjusting

Use Touch Panel Trackball to position the M cursor over the required area of the B-Mode image.

1. Use the trackball (assigned function: Pos) to position the M cursor over the required area of the image.

Performing an Exam

2. Press the top Trackball key to allow free rotation of the solid arrow line throughout the 2D image (trackball assigned function: Angle).

NOTE

Rotate the Touch Panel control to angle the M cursor at a given angle.

Benefits

Color Flow Mode and Color M-Mode are Doppler Modes intended to add color-coded qualitative information concerning the relative velocity and direction of fluid motion within the B-Mode or M-Mode image.

Bioeffects

Changing the Packet Size, Scale, and ROI size may change the TI and/or MI. Observe the output display for possible effects.

Curved Anatomical M-Mode (Camm)

Curved Anatomical M-Mode (Camm) displays a distance/time plot from a free-drawn cursor line. Camm is available in B, CF and TVI.

1. Select **Curved AMM** on the Touch Panel.
2. Use the Trackball to position the start point of the time motion curve in the B-Mode image.
3. Press **Set** to fix the start point.
4. Use the Trackball to position the next point.
The time motion curve is drawn by the green line.
5. Press **Set** to fix the point.
6. Repeat step 4 and 5 to draw a complete time motion curve.

NOTE

The time motion curve can be edited by following the curve back to the desired point and redrawn as desired. Following the curve back to the starting point will delete the time motion curve.

7. Press **Set** twice to complete.

NOTE

Move the cursor to the desired anchor point and press **Set**. Move the point to the desired position and press **Set**.

8. The arrow cursor appears on the M-Mode image and the red bar appears on the time motion curve.

The red bar indicates the position of the time motion curve relative to the arrow cursor on the Camm image. They move relative to one another.

NOTE

Press **Set** to clear a cursor line.

NOTE

Curved Anatomical M-Mode can also be used with previously acquired digitally stored B-Mode images.

NOTE

CAMM is not available on Linear probes.

Optimizing Color Flow Mode

Color Flow Mode is a Doppler Mode intended to add color-coded qualitative information concerning the relative velocity and direction of fluid motion within the B-Mode image.

Typical exam protocol

A typical examination using Color Flow Mode.

1. Follow the same procedure as described under B-Mode to locate the anatomical area of interest.
2. After optimizing the B-Mode image, add Color Flow.

NOTE

Use all noise reduction controls with care. Excessive application may obscure low level diagnostic information.

3. Move the color flow area of interest as close to the center of the image as possible.
4. Optimize the color flow parameters so that a high frame rate can be achieved and appropriate flow velocities are visualized.
5. Press Freeze to hold the image in memory.
6. Record color flow images as necessary.
7. If more definitive information is needed about flow, utilize the procedures described under Doppler Mode.
8. To exit Color Flow, press CF-Mode or B-Mode.

NOTE

Most parameters are user presettable by probe and application in the preset menu (*Utility > Imaging > CF*).

Color Flow Mode Controls

Color Flow Mode and Color M-Mode are Doppler Modes intended to add color-coded qualitative information concerning the relative velocity and direction of fluid motion within the B-Mode or M-Mode image.

Table 2-3 Flow Mode Controls

Control	Bioeffect	Description/Benefit
Flow Selection	Yes	In the Lower Extremity Vein (LEV) and Abdominal applications, you can quickly select the flow state via a shortcut on the Color Flow Mode menu.

Control	Bioeffect	Description/Benefit
Gain	No	Gain amplifies the overall strength of echoes processed in the Color Flow window or spectral Doppler timeline.
Scale (Velocity Scale)	Yes	Increases/decreases the Scale on the color bar.
Wall Filter	No	Filters out low flow velocity signals. It helps get rid of motion artifacts caused from breathing and other patient motion.
Size/Position	Yes	Adjust size and position of the color window.
Invert (Color Invert)	No	Lets you view blood flow from a different perspective, e.g., red away (negative velocities) and blue toward (positive velocities). You can invert a real-time or frozen image. NOTE Invert reverses the color map, NOT the color PRF.
Baseline	No	Changes the Color Flow or Doppler spectrum baseline to accommodate higher velocity blood flow. Minimizes aliasing by displaying a greater range of forward flow with respect to reverse flow, or vice versa. Baseline adjusts the alias point. The default baseline is at the midpoint of the color display and at the midpoint of the color bar reference display.
Angle Steer	Yes	You can slant the ROI of the Color Flow linear image left or right to get more information without moving the probe. The Angle Steer function only applies to linear probes.
Accumulation	No	Time Accumulation of Color Flow.
Color Flow Line Density	Yes	Optimizes the Color Flow frame rate or spatial resolution for the best possible color image.
Map	No	Allows you to select a specific color map. After you have made your selection, the color bar displays the resultant map.
Map Compress	No	Change the gradation of color map.
Threshold	No	Limits color flow overlay to low level echoes inside vessel walls. Helps minimize color 'bleeding' outside vessel walls.
Frame Average	No	Averages color frames.
Transparency Map	No	Brings out the tissue behind the color map.
Spatial Filter	No	Smooths out the color and reduces over pixelated appearance.
Flash Suppression	No	Activates/deactivates Flash Suppression, a motion artifact elimination process.
Packet Size	Yes	Controls the number of samples gathered for a single color flow vector.

Control	Bioeffect	Description/Benefit
Sample Volume	Yes	Adjusts the size of the color flow doppler transmit wave (or pulse) and size (or length). Lower setting gives better flow resolution and a higher setting increases sensitivity.
CF/PDI Auto Sample Volume	Yes	Set the default value at Utility > Imaging > CF Mode. The transmit frequency and focus position will be changed automatically with depth.
CF/PDI Focus Depth	Yes	
CF/PDI Frequency	Yes	
CF/PDI Auto Frequency	Yes	
CF/PDI Center Depth	Yes	
PDI	No	Power Doppler Imaging (PDI) is a color flow mapping technique used to map the strength of the Doppler signal coming from the tissue rather than the frequency shift of the signal. Using this technique, the ultrasound system plots color flow based on the number of reflectors that are moving, regardless of their velocity. PDI does not map velocity, therefore it is not subject to aliasing.
Directional PDI	No	You can select the DP7 Directional Power Doppler maps while in PDI. Directional Maps. Directional Power Doppler Maps from PDI. NOTE If you store a PDI image and recall it, you can still switch to the Directional Power Doppler map and vice versa. However, an image stored as non-directional then switched to directional just adds direction to a non-directional map and vice versa.
TVI (Option)	Yes	Tissue Velocity Imaging (TVI) calculates and color-codes the velocities in tissue. The tissue velocity information is acquired by sampling of tissue Doppler velocity values at discrete points. The information is stored in a combined format with gray scale imaging during one or several cardiac cycles with high temporal resolution.
TVD (Option)	Yes	Tissue Velocity Doppler (TVD): basing on TVI mode, activate a sample volume of PW ventricular wall to get the spectral information of the sample section.
TVM (Option)	Yes	Tissue Velocity M Mode image (TVM): Active M mode while in TVI. To display myocardium motion velocity and direction.
QAnalysis (Option)	Yes	QAnalysis is available for the image loop acquired in the following modes: TVI, CF and PDI. All of the QAnalysis modes operate similarly, with some variation.

Optimizing M Color Flow

M Color Flow is used for fetal cardiac applications. Color Flow overlays color on the M-Mode image using velocity and variance color maps. The Color Flow wedge overlays the B-Mode image and M-Mode timeline.

The Color Flow maps available in M-Mode are the same as in Color Flow Mode. The size and position of the Color Flow window in B-Mode determines the size and position of the Color Flow window in M-Mode.

All M-Mode measurements are available with M Color Flow active: depth, distance along a straight line, % stenosis, volume, trace, circumference, enclosed area, distance, time, slope, and heart rate.

M Color Flow mode can be activated on the sector probes only.

Activating

To activate M Color Flow Mode, press **M** (M-Mode), then press **CF** (Color Flow)- or - press **CF**, then press **M**.

CM tab is displayed.

Benefits

Color Flow Mode and Color M-Mode are Doppler Modes intended to add color-coded qualitative information concerning the relative velocity and direction of fluid motion within the B-Mode or M-Mode image.

Bioeffects

Changing the Sweep Speed, Packet Size, Frame Rate/Resolution, Zoom, PRF, and ROI size may change the TI and/or MI. Observe the output display for possible effects.

Optimizing Doppler Mode

Doppler is intended to provide measurement data concerning the velocity of moving tissues and fluids. PW Doppler lets you examine blood flow data selectively from a small region called the sample volume.

Typical Use - PW Doppler

In Pulsed Wave Doppler (PW) Mode, energy is transmitted from the ultrasound probe into the patient, as in B-Mode. However, the received echoes are processed to extract the difference in frequency between the transmitted and received signals. Differences in frequencies can be caused by moving objects in the path of the ultrasound signal, such as moving blood cells. The resultant signals are presented audibly through the system speakers and graphically on the system display. The X axis of the graph represents time while the Y axis represents the shift in frequency. The Y axis can also be calibrated to represent velocity in either a forward or reverse direction.

PW Doppler is typically used for displaying the speed, direction, and spectral content of blood flow at selected anatomical sites. PW Doppler operates in two different modes: conventional PW and High Pulse Repetition Frequency (HPRF).

PW Doppler can be combined with B-Mode for rapidly selecting the anatomical site for PW Doppler examination. The site where PW Doppler data is derived appears graphically on the B-Mode image (Sample Volume Gate). The sample volume gate can be moved anywhere within the B-Mode image.

Typical exam protocol

A typical examination using PW Doppler Mode might proceed as follows:

1. Locate the anatomy to be examined. Get a good B Mode image. Press **CF** to help locate the vessel you wish to examine.
2. Press **PW**. The PW Doppler spectrum appears and the system operates in combined B mode, color Doppler and PW Doppler modes. Adjust **Volume** to adjust Doppler audio. The Doppler signal is heard through the speakers.
To activate Continuous Wave (CW), press **CW**.
3. Position the sample volume cursor by moving the **Trackball** left and right. Position the sample volume gate by moving the **Trackball** up and down. Size the gate by clicking **SV Length**.
4. Optimize the PW Doppler spectrum, as necessary. Refer to the Doppler Optimization section of this chapter for more information.
5. Press **Set** to toggle between real time B-Mode with Doppler Mode (with audio).
6. Sample along the whole length of the vessel. Make sure that the probe is parallel to flow. Listen, then look, when positioning the sample volume cursor.
7. Press **Freeze** to hold the trace in memory and stop imaging. Activate CINE Timeline, as necessary.
8. Perform measurements and calculations, as necessary. Refer to the Measurements and Calculations chapter for more information.
9. Record results by pressing the appropriate print key, depending on the setup of your recording devices.
10. Press Freeze to resume imaging.
11. Repeat the above procedure until all relevant flow sites have been examined.

Doppler Mode Display Explanations

Table 2-4 Doppler Mode Display Explanations

Doppler Display	Description, Format, Values
PRF	Pulse repetition frequency, displayed as PRF in kHz.
Wall Filter	Wall filter size, displayed as WF in Hz.
Doppler Gain	Displays as GN in decibels (dB).
Sample Volume Depth	Displays (in cm) when Doppler cursor is present.
Doppler Angle (AC ##)	Indicates angle in degrees between the Doppler mode cursor and the angle correction indicator. Displays when Doppler cursor is present. The Doppler Angle displays in red when the angle exceeds 60°. Velocities obtained when the angle is greater than 80° are displayed as asterisks (**).

Doppler Display	Description, Format, Values
Spectral Invert	INVERT appears when the spectral trace is inverted and the plus/minus signs (+/-) are reversed.
HPRF	HPRF mode is used when detected velocities exceed the processing capabilities of the currently selected PW Doppler scale or when the selected anatomical site is too deep for the selected PW Doppler scale.
Time Scale	Each selection represents a different sweep time.
Angle Correct	Indicates flow direction.
Sample Volume Gate	Indicates sample volume box. Each probe defaults to a specific range gate.
Doppler Velocity Scale	Flow direction has a positive and negative indicator, noted in centimeters per second (cm/sec.). When the velocity scale is less than 10 cm/s, it is displayed to the first decimal point (4.6 rather than 5 cm/s). The Doppler velocity scale adjust as you adjust the PRF.

Doppler Mode Controls
Table 2-5 Doppler Mode Controls

Control	Bioeffect	Description/Benefit
Auto Spectral Optimize [ASO] (Auto)	Yes	Auto in Doppler Mode optimizes the spectral data. Auto adjusts the Velocity Scale/PRF (on live images only), baseline shift and invert (if preset). The benefit of Auto can be found in reduced optimization time and a more consistent and accurate optimization process.
Set	Yes	Toggles between simultaneous and update presentation while viewing the timeline.
Doppler sample volume gate position	Yes	Moves the sample volume gate on the B-Mode's Doppler Mode cursor. The gate is positioned over a specific position within the vessel.
Doppler sample volume length	Yes	Sizes the sample volume gate.
Scale (Velocity Scale)	Yes	Adjusts the velocity scale to accommodate faster/slower blood flow velocities. Velocity scale determines pulse repetition frequency. If the sample volume gate range exceeds single gate Scale capability, the system automatically switches to high PRF mode. Multiple gates appear, and HPRF is indicated on the display.

Control	Bioeffect	Description/Benefit
Angle Correct	No	Estimates the flow velocity in a direction at an angle to the Doppler vector by computing the angle between the Doppler vector and the flow to be measured. NOTE When the Doppler Mode Cursor and angle correct indicator are aligned (the angle is 0), you cannot see the angle correct indicator.
Quick Angle	No	Quickly adjusts the angle by 60 degrees.
Wall Filter	No	Insulates the Doppler signal from excessive noise caused from vessel movement.
Baseline	No	Adjusts the baseline to accommodate faster or slower blood flows to eliminate aliasing.
Mode Cursor	No	Displays the Doppler Mode cursor on the B-Mode image.
Steer and Fine Steer	Yes	You can slant the ROI of the Color Flow linear image left or right to get more information without moving the probe. The angle steer function only applies to linear probes.
Volume	No	Controls audio output.
Sample Volume	Yes	Adjusts the size of the color flow doppler transmit wave (or pulse) and size (or length). Lower setting gives better flow resolution and a higher setting increases sensitivity.
Invert	No	Vertically inverts the spectral trace without affecting the baseline position.
Compression	No	Dynamic range controls how echo intensities are converted to shades of gray, thereby increasing the range of contrast you can adjust.
Trace Method	No	Traces the average mean and peak velocities in real-time or frozen images.
Cycles to Average	No	The average value over a number of cycles (from 1-5).
Trace Sensitivity	No	Adjust the trace to follow the waveform for signal strength.
Trace Direction	No	Specifies trace direction.
Display Format	No	Changes the horizontal/vertical layout between B-Mode and Doppler-Mode, or timeline only.
Modify Auto Calcs	No	Activates the menu to select which calculations are automatically calculated.
Auto Calcs	No	Activates the calculation automatically which you select in the Modify Auto Calculation when the system is in a state of freeze or live.

Control	Bioeffect	Description/Benefit
Simultaneous (Duplex/Triplex)	Yes	When you select Simultaneous, everything is live. Simultaneous includes duplex and triplex. For example, it is duplex if both B-Mode and PW Doppler Modes are active; it is triplex if B-Mode, PW Doppler Mode, and CF Modes are active. If Simultaneous is not selected, use Set to toggle between modes.
Continuous Wave Doppler (CWD) (option)	Yes	Allows examination of blood flow data all along the Doppler Mode cursor rather than from any specific depth. Gather samples along the entire Doppler beam for rapid scanning of the heart. Range gated CW allows information to be gathered at higher velocities.
Frequency	Yes	In CWD and PWD mode, touch the frequency button on the touch screen to adjust PWD and CWD mode frequency.

Other Controls

Zoom

To zoom an image, rotate / press **Zoom** knob. A reference image appears in the lower left corner of the monitor.

To exit zoom, enter **B** mode or press **Zoom** knob until the reference zoom image is removed. Reference Image is the small un-zoomed image displayed next to the zoomed image.

Read Zoom

System will enter **Read Zoom** by default after press / rotate **Zoom** knob.

Read Zoom magnifies the display of the data without making any changes to the ultrasound image data that is acquired.

Available in a live, frozen, cine or recalled raw data image.

Write Zoom

To activate Write Zoom, press the trackball key to enter **Write Zoom**.

Zooming an image changes the frame rate which tends to change thermal indices. The position of the focal zones may also change which may cause the peak intensity to occur at a different location in the acoustic field. As a result, the MI (TI) may change.



ACOUSTIC OUTPUT HAZARD

Observe the output display for possible effects.

With Write Zoom, the Ultrasound line density and/or sampling frequency increases, giving a better resolution.

Write Zoom is available only in pre-processing.

You can preset the write zoom window size (height and width) on **Utility > Imaging > B-Mode**.

NOTE

The difference between Read Zoom and Write Zoom can be described in relation to photography. With a photograph, Read Zoom manipulates the negative and enlarges the picture; whereas Write Zoom uses a telephoto lens to bring the image closer before taking the picture.

NOTE

First use the Read Zoom (turn knob) to get to the area of interest, then use Write Zoom (press knob).


V-Zoom

V-Zoom is intended to zoom the full screen for selected ROI. It works for both **Read Zoom** and **Write Zoom**.

To activate V-Zoom, follow below steps:

1. Select **Enable vZoom** through **Utility > System > System Imaging**.

Figure 2-6 Enable vZoom in Utility




2. Press/rotate the **Zoom** knob during scanning, an ROI frame displayed on the screen.

Performing an Exam


3. Adjust position and size by Pos/Size key.

Figure 2-7 Adjust position and size



4. Press the **Zoom** knob again to display the magnified ROI.

Figure 2-8 V-Zoom example



To exit V-Zoom, press **B** mode or the **Zoom** knob.

Split screen

To activate a dual split screen, press Dual. To activate a quad display, press Quad.

To switch between active images, press Dual/Quad.

To deactivate, press Quad until the screen changes.

NOTE

To put a copy of the image on the opposite side when entering dual split screen, use the ***When Entering Dual Image...*** preset via ***Utility > Application > Settings***.

Freezing an image

To freeze an image, press **Freeze** key. The **Freeze** key backlight turns green. Image is on freeze mode (Frozen image).

If you are in a mixed mode, both images freeze immediately. Deactivating Freeze restarts both modes.

To reactivate the image, press **Freeze** key again. The **Freeze** key backlight becomes blue (unfreeze). Image is on real time scan mode (Live image).

NOTE

Deactivating Freeze erases all measurements and calculations from the display (but not from the report).

Use the Trackball to start CINE after pressing Freeze.

Activating CINE Mode

Follow below steps to activate CINE Mode.

1. Press **Freeze** key.
2. Move the Trackball to activate CINE.

Body Patterns

Press the Body Pattern key on the control panel to enter into body pattern mode.

Select the desired body pattern and the selected body pattern is displayed on the monitor.

- Select the **Move Pattern** control to reposition the body pattern with the **Trackball** and **Set** controls.
- Move the body pattern to the desired location.

NOTE

Body Pattern Position is reset to factory default when patient is changed (e.g. End Exam).

- A probe mark is associated with the body patterns and illustrates the probe position on the body pattern. This marker can be placed with the **Trackball** and rotated with the **Body Pattern** button.
- The probe mark type is selectable by rotating the **Probe Type** control. There are different choices available with one being a blank selection.
- To select the active side in dual B-Mode, use the **Active Side** rotary control on the Primary Menu.
- To clear the body pattern, press the **Body Pattern** control to activate body patterns and then press the **Clear** key.
- Press **Set** on the keyboard to exit without erasing the body pattern.

Annotating an Image

Pressing the Comment key or any keys on the alphanumeric keyboard initiates the comment mode. This assigns the trackball function to controlling the cursor and displays the comment library on the menu area.

In comment mode, text can be added by using the comment library or by typing from the alphanumeric keyboard.

After activating the comment mode, a vertical bar type cursor appears on the screen. Use the Trackball to move the cursor.

To delete comments by character, press the Backspace key.

To delete all comments and arrow marks, press the Clear key after entering the comment mode.

To move by words or by text group, press the Tab key.

Arrow pointers can be used by activating the Arrow key on the keyboard. Comment key can be activated either by pressing the Annotate key or Arrow key on the control panel. When the pointer comes up, it is a GREEN color, indicating it is active and can be moved.

NOTE

Please notify that there is limitation on the length of annotations.

Measurement and Analysis


Measurements and calculations derived from ultrasound images are intended to supplement other clinical procedures available to the attending physician. The accuracy of measurements is not only determined by system accuracy, but also by the use of proper medical protocols by the user. When appropriate, be sure to note any protocols associated with a particular measurement or calculation. Formulas and databases used within the system software that are associated with specific investigators are so noted. Be sure to refer to the original article describing the investigator's recommended clinical procedures.

NOTE

The system provides calculations (e.g. estimated fetal weight) and charts based on published scientific literature. The selection of the appropriate chart and clinical interpretation of calculations and charts is the sole responsibility of the user. The authorized user should consider proper indications for the use of a calculation or chart as described in the scientific literature. The diagnosis, decision for further examination, and medical treatment must be performed by qualified personnel following good clinical practice.

Location of Measurement Controls

Figure 2-9 Locating Measurement Controls - example



1. **Cursor** - Pointer Key. Select to display a pointer on the monitor.
2. **Clear** - During a measurement sequence, erases the measuring caliper and measurement data from the display. When not performing a measurement sequence, clears all calipers and measurements from the display.
3. **Measure** - Activates a measurement caliper and the calculation package associated with the currently selected preset.

4. **Body Pattern** - Body patterns can be added to measurements on the image.
5. **Zoom/Ellipse**
 - To Zoom an image, rotate **Zoom** clockwise or press **Zoom** knob. For more details, refer to *Zoom* on page 90.
 - Press **Ellipse** knob to activate ellipse function. After the first caliper for a distance measurement has been set and the second caliper positioned, Ellipse activates the area/ellipse measurement function. To increase the size, adjust the **Ellipse** control in a clockwise direction. To decrease the size, adjust the **Ellipse** control in a counterclockwise direction.
6. **Comment** - Comments can be added to measurements on the image.
7. **Trackball** - Moves the measurement calipers, selects the measurement on the Summary Window. Trackball also selects items on the Touch Panel with the Pointer and Set keys.
8. **Trackball keys** - The functionality of these keys changes (for example, Set, Change Measure, etc.) depending on the mode or action. Current functionality is displayed on the lower right corner of the monitor.

B-Mode Measurements

The following measurements can be made in B-Mode.

- Distance
- Circumference
- Circumference and Area
 - Ellipse Method
 - Trace Method
 - Spline Method
 - Intensity (Echo level) Method

NOTE

Prior to performing B-Mode measurements, scan the patient and then press **Freeze**.



WARNING

DO NOT perform a depth measurement using 4D probes.

Distance Measurement

To make a distance measurement:

1. Press **Measure** once; an active caliper displays.
2. To position the active caliper at the start point, move the **Trackball**.
3. To fix the start point, press **Set**.

The system fixes the first caliper and displays a second active caliper.
4. To position the second active caliper at the end point, move the **Trackball**.

A dotted line connects the measurement points, if preset accordingly.

5. To complete the measurement, press **Set**.

The system displays the distance value in the Results Window.



HINTS

- **Before** you complete a measurement:
 - To toggle between active calipers, press the top Trackball key.
 - To erase the second caliper and the current data measured and start the measurement again, press **Clear** once.
- **After** you complete the measurement:
 - To rotate through and activate previously fixed calipers, adjust **Cursor Select**.
 - To erase all data that has been measured to this point, but not data entered onto worksheets, press **Clear**.

Circumference Measurement (Open Trace)

To trace the circumference of a portion of the anatomy and calculate its length:

NOTE

Set Open Trace to the Touch Panel in **Utility > Measure** before perform the measurement.

1. Press **Measure**.
2. Select **Open Trace** from the Touch Panel.
3. Position the caliper at the start point.
4. To fix the trace start point, press **Set**. The caliper changes to an active caliper.
5. Move the **Trackball** to trace the measurement area. A dotted line shows the traced area.
6. To complete the measurement, press **Set**. The system displays the circumference in the Results Window.

Circumference and Area (Ellipse) Measurement

You can use an ellipse to measure circumference and area. To measure with an ellipse:

1. Press **Measure** once; an active caliper displays.
2. To position the active caliper, move the **Trackball**.
3. To fix the start point, press **Set**. The system fixes the first caliper and displays a second active caliper.
4. To position the second caliper, move the **Trackball**.
5. Rotate the **Ellipse** control; an ellipse with an initial circle shape displays.

6. To position the ellipse and to size the measured axes (move the calipers), move the **Trackball**.
7. To increase the size, rotate the **Ellipse** control in a clockwise direction. To decrease the size, rotate the **Ellipse** control in a counterclockwise direction.
8. To toggle between active calipers, press the top **Trackball key**.
9. To complete the measurement, press **Set**. The system displays the circumference and area in the Results Window.



HINTS

Before completing the ellipse measurement:

- To erase the ellipse and the current data measured, press **Clear** once. The original caliper is displayed to restart the measurement.
- To exit the measurement function without completing the measurement, press **Clear** a second time.

Circumference and Area (Trace) Measurement

To trace the circumference of a portion of the anatomy and calculate its area:

1. Press **Measure**.
2. Press the top **Trackball key** to select Trace; a caliper displays.
3. To position the caliper at the start point, move the **Trackball**.
4. To fix the trace start point, press **Set**. The caliper changes to an active caliper.
5. To trace the measurement area, move the **Trackball** around the anatomy. A dotted line shows the traced area.
6. To complete the measurement, press **Set**. The system displays the circumference and the area in the Results Window.



HINTS

Before completing the trace measurement:

- To erase the line (bit by bit) back from its current point, move the **Trackball** or adjust the **Ellipse** control counterclockwise.
- To erase the dotted line but not the caliper, press **Clear** once.
- To clear the caliper and the current data measured, press **Clear** twice.

Circumference and Area (Spline Trace) Measurement

To trace the circumference of a portion of the anatomy and calculate its area:

NOTE

Spline trace is not available through the factory default. The system defaults to trace. To enable spline trace, modify the Measure Key Sequence preset found in Utility > Measure > Advanced preset menu.

1. Press **Measure**.

2. Press the top **Trackball key** to select Spline Trace; a caliper displays.
3. To position the first caliper at the start point, move the **Trackball**.
4. To fix the trace start point, press **Set**. The first caliper turns yellow. The second caliper appears at the same position as the first caliper and is green.

NOTE

When pressing the **Clear** key once, the second caliper disappears and the first caliper is activated.

If **Clear** is pressed again, the first caliper disappears and the Spline trace is cancelled.

5. To position the second caliper, move the **Trackball** and press **Set**. The third caliper appears at the same position.

NOTE

The **Clear** key functionality is the same as noted in the previous step.

The spline trace requires at least three points to draw the trace. Continue setting the points of the trace until the desired points are set.

6. Press **Set** again after the last caliper is fixed to finalize the spline trace. All points are removed from the line and the spline trace turns yellow.

NOTE

Pressing **Set** twice finishes the trace measurement.

If **Clear** is pressed twice when more than 3 points exist on the trace, all points are removed and the first caliper again displays.

Edit the Spline Trace


1. Select **Cursor Select**. The spline trace changes to green and all points appear on the trace as yellow.

A pick-caliper appears on the center of the image and the message "Edit spline trace" displays at the bottom of the screen.

NOTE

The pick-caliper is used to select and move the trace points.

Figure 2-10 Edit spline trace



Select **Cursor Select** again. The trace is deactivated (changes to yellow) and all points, including the pick-caliper, are removed.

If the previous/next fixed caliper exists on the image, it is activated.

NOTE

Pressing **Clear** at this time removes all points and the trace graphic.

2. Move the pick-caliper to the desired point and press **Set**. The point is activated and turns green.

Performing an Exam

3. Move the point to the desired position and press **Set**. The point is fixed and turns yellow. The pick-caliper appears on the center of the image.

NOTE

The spline trace is updated while system is running.

NOTE

To remove a point, press **Clear** while moving the point. The trace turns green and the remaining points continue to be shown as yellow. If there are less than three points, the spline trace is removed.

4. Press **Set** again. All points are removed from the trace and the trace is shown as yellow.

Intensity (Echo Level) Measurement

To make an echo level measurement:

1. Press **Measure**.
2. Press the top Trackball key to select Intensity. A caliper displays.
3. To position the caliper at the start point, move the **Trackball**.
4. To fix the trace start point, press **Set**. The caliper changes to an active caliper.
5. To trace the measurement area, move the **Trackball** around the anatomy. A dotted line shows the traced area.
6. To complete the measurement, press **Set**. The system displays the echo level, as EL __ dB, in the Results Window.

NOTE

The echo level measurement is only available on a frozen image, not on a B-paused image.

NOTE

Echo Level is not available through the factory default. To enable echo level, modify the Measure Key Sequence preset, found in the Utility > Measure > Advanced preset.

M-Mode Measurements

Basic measurements that can be taken in the M-Mode portion of the display are:

- Tissue Depth (Distance)
- Time Interval
- Time Interval and Velocity



WARNING

DO NOT perform a depth measurement using 4D probes.

Prior to performing M-Mode measurements:

1. In the B-Mode part of the display, scan the anatomy you want to measure.
2. Go to the M-Mode part of the display.
3. Press **Freeze**.

Tissue Depth

Tissue depth measurement in M-Mode functions the same as distance measurement in B-Mode. It measures the vertical distance between calipers.

1. Press **Measure** once; an active caliper with a vertical and horizontal dotted line displays.
2. To position the active caliper at the most anterior point you want to measure, move the **Trackball**.
3. To fix the start point, press **Set**.
The system fixes the first caliper and displays a second active caliper.
4. To position the second caliper at the most posterior point you want to measure, move the **Trackball**.
5. To complete the measurement, press **Set**.

The system displays the vertical distance between the two points in the Results Window.

Time Interval

To measure a horizontal time interval and velocity:

1. Press **Measure**. Press the top Trackball key to select Time; an active caliper with vertical and horizontal dotted lines displays.
2. To position the caliper at the start point, move the **Trackball**.
3. To fix the first caliper, press **Set**. The system fixes the first caliper and displays a second active caliper.
4. To position the second caliper at the end point, move the **Trackball**.
5. To complete the measurement, press **Set**. The system displays the time interval between the two calipers in the Results Window.

Slope (Time Interval and Velocity)

To measure time and velocity between two points:

1. Press **Measure**. Press the top Trackball key to select Slope; an active caliper with vertical and horizontal dotted lines displays.
2. To position the active caliper at the start point, move the **Trackball**.
3. To fix the start point, press **Set**.
The system fixes the first caliper and displays a second active caliper.
4. To position the second caliper at the end point, move the **Trackball**.
5. To complete the measurement, press **Set**.

The system displays time(s) and slope between the two points in the Results Window.

Doppler Mode Measurements

Five basic measurements can be made in Doppler Mode.

- Velocity
- TAMAX and TAMEAN (Manual or Auto Trace)
- Two Velocities with the Time Interval and Acceleration between them
- Time Interval
- Volume Flow

Prior to performing Doppler Mode Measurements:

1. In the B-Mode part of the display, scan the anatomy you want to measure.
2. Go to the Doppler Mode part of the display.
3. Press **Freeze**.

Velocity

To measure velocity:

1. Press **Measure**; an active caliper with a vertical dotted line displays.
2. To position the caliper at the desired measurement point, move the **Trackball**.
3. To complete the measurement, press **Set**. The system displays the velocity measurement in the Results Window.

TAMAX and TAMEAN Manual Trace

The value measured depends upon the Vol Flow Method preset. The two selections available are: Peak (TAMAX) and Mean (TAMEAN).

To do a manual trace of TAMAX or TAMEAN:

1. Press **Measure**. Press the top Trackball key to select Trace; a caliper displays. Select Manual on the Touch Panel.
2. To position the caliper at the trace start point, move the **Trackball**.
3. To fix the start point, press **Set**.
4. To trace the velocity spectrum boundary, move the **Trackball**.

NOTE

To edit the trace line, move the **Trackball**.

5. To complete the measurement, press **Set**. The system displays the measurement values in the Results Window.

Auto Trace

The value measured depends upon the Vol Flow Method preset. The two selections available are: Peak (TAMAX) and Mean (TAMEAN).

To auto trace TAMAX:

1. Press **Measure**. Press the top Trackball key to select Trace; an active caliper with a vertical dotted line displays. Select Auto on the Touch Panel.
2. To position the caliper at the trace start point in the Doppler spectrum, move the **Trackball**.
3. To fix the start point, press **Set**.
4. To position the vertical caliper at the end point, move the **Trackball**.
5. To complete the measurement, press **Set**. The system automatically fixes both calipers and traces the maximum value between the two points. The system displays this value in the Results Window.

NOTE

When you set the Auto Trace for Both (above and below), the system picks up the maximum power of the signal, NOT the maximum velocity. If the maximum velocity is not the maximum power, the system may not trace accurately. If you want to use maximum velocity, select either Above or Below.

Edit Trace

Auto Trace can be edited after taking an Auto Trace measurement.

1. After taking an Auto Trace measurement, select the measurement result on the result window. The Edit Trace (Edit Peak or Edit Mean) menu window appears.

NOTE

If the system cannot take the trace data correctly from the image, Edit Trace does not work.

2. Select Edit Trace. The first caliper (manual trace caliper) appears on the center of the image. Use the **Trackball** to move the caliper on the trace line to the start point.

NOTE

To cancel Edit Trace at this time, press **Clear**, **Scan**, or **Freeze**.

3. Press **Set** to fix the first caliper. The second caliper appears. Edit the trace manually using the second caliper.

The Ellipse control is used to edit the trace.

NOTE

When pressing the **Clear** key once at this time, the second caliper disappears and the first caliper appears in the center of the image.

NOTE

If you press **Scan** or **Freeze** at this time, the caliper is automatically fixed and the result window updates.

4. Press **Set** to fix the second caliper. The trace and the result window update. The trace data (TAMAX and TAMEAN) are updated, though the other points (e.g. PS, ED) are not updated by trace. The points can be edited with **Cursor Select**.

NOTE

While in Edit Trace, Cursor Select is disabled.

5. Repeat Edit Trace as needed.

Slope (Velocity, Time Interval and Acceleration)

To measure two velocity values, the time interval (ms), and acceleration (m/s^2):

1. Press **Measure**. Press the top Trackball key to select Slope; an active caliper with vertical and horizontal dotted lines displays.
2. To position the caliper at the start point, move the **Trackball**.
3. To fix the start point, press **Set**. The system fixes the first caliper and displays a second active caliper.
4. To position the second caliper at the end point, move the **Trackball**.
5. To complete the measurement, press **Set**. The system displays the two peak end point velocities, the time interval, and the acceleration in the Results Window.

Time interval

To measure a horizontal time interval:

1. Press **Measure**. Press the top Trackball key to select Time; an active caliper with vertical and horizontal dotted lines displays.
2. To position the active caliper at the start point, move the **Trackball**.
3. To fix the start point, press **Set**. The system fixes the first caliper and displays a second active caliper.
4. To position the second caliper at the end point, move the **Trackball**.
5. To complete the measurement, press **Set**. The system displays the time interval between the two calipers in the Results Window.

Volume

The volume calculation can be made from any of the following measurements:

- One distance
- Two distances
- Three distances
- One ellipse
- One distance and one ellipse

For details on how to make a distance measurement, refer to *Distance Measurement* on page 96.

For details on how to make an ellipse measurement, see *Circumference and Area (Ellipse) Measurement* on page 97.

NOTE

IMPORTANT!! If you want to make a volume calculation using one or two distances, you must select **Volume** BEFORE you make the measurements.

NOTE

If you select Fix Caliper by Print Key on the **Utility > System > System Measure**, the print key does not function like the **Set** key, but instead ends the measurement sequence and initiates the volume calculation based on the number of measurements taken so far.

To make a volume calculation using one or two distances:

1. Select **Volume**.
2. Make one or two distance measurements.
3. Select **Volume**.
The system displays the distances and the volume in the Results Window.

NOTE

Use the **Clear** key to erase the green caliper.

To make a volume calculation using three distances:

1. Make three distance measurements.

NOTE

Three distances can be done in the dual format mode (side by side images). One measurement is usually made in the sagittal plane and two measurements in the axial plane. To use the dual format mode, press the **L** or **R** key on front panel.

2. Select **Volume**.
The system displays the distances and the volume in the Results Window.

To make a volume calculation using one ellipse:

1. Make one ellipse measurement.
2. Select **Volume**.
The system displays the ellipse measurement and the volume in the Results Window.

To make a volume calculation using one ellipse and one distance:

1. Make one distance measurement and one ellipse measurement.
2. Select **Volume**.
The system displays the distance and ellipse measurement and the volume in the Results Window.



HINTS

- Volumes are most accurate when measurements are taken in the sagittal and axial scan planes.
- To display sagittal and axial plane images simultaneously, use the side-by-side dual format option.

NOTE


If you change the parameters or category during the volume measurement, please follow the procedure below before you restart the measurement.

1. Check the number of each measurement in the summary window.
2. If the numbers are not all the same, it shows that you have the calculation which is not completed. Open the Worksheet and clear that calculation.

Table 2-6 Volume Calculations

Calc Name	Input Measurements
Volume (spherical)	One distance
Volume (prolate spheroidal)	Two distances, $d1 > d2$
Volume (spheroidal)	Three distances
Volume (prolate spheroidal)	One ellipse: (d1 major axis, d2 minor axis)
Volume (spheroidal)	One distance d1, and one ellipse (d2 major axis, d3 minor axis)

Figure 2-11 Volume Calculation Examples




Calculation formulas are available in the Advanced Reference Manual.

Post-Assignment for General Volume

You can input a unique name for the general volume measurement. You can group the general volume measurements for each application.


1. Complete the volume measurement.
2. Move the caliper to the measurement result box (with green frame) and select **Set**.
3. The volume name menu appears. Select **Name Volume**.

Figure 2-12 Volume Name Menu



4. The dialog box displays. Enter a new name or choose the existing name.

Figure 2-13 Volume Name Dialog box



NOTE


The factory default volume name cannot be changed (for example, Renal Volume).

Application M&A

Abdomen

Abdominal measurements offer a few different types of measurement studies:

Figure 2-14 Abdomen Exam Category



- Generic-Common to all applications.
- Abdomen
- Renal
- Bypass Graft
- Aorta Iliac
- Vascular Renal
- Mesenteric
- Abdomen Vein
- Bladder

To change a study:

1. Press **Measure**.
2. Select **Abdomen** in **Exam Category**.
3. To choose another study, select the desired study.

Small Parts

B-Mode Measurements

The Small Parts exam category includes the following folders:

- Generic - Common to all applications.
- Thyroid
- Scrotal

- Breast

OB

Prior to an ultrasound examination, the patient should be informed of the clinical indication, specific benefits, potential risks, and alternatives, if any. In addition, if the patient requests information about the exposure time and intensity, it should be provided. Patient access to educational materials regarding ultrasound is strongly encouraged to supplement the information communicated directly to the patient. Furthermore, these examinations should be conducted in a manner and take place in a setting which assures patient dignity and privacy.

- Prior material knowledge and approval of the presence of nonessential personnel with the number of such personnel kept to a minimum.
- An intent to share with the parents per the physician's judgement, either during the examination or shortly thereafter, the information derived.
- An offer of choice about viewing the fetus.
- An offer of choice about learning the sex of the fetus, if such information becomes available.

Ultrasound examinations performed solely to satisfy the family's desire to know the fetal sex, to view the fetus, or to obtain a picture of the fetus should be discouraged.

Acoustic Output Considerations



DANGER

The ultrasound system is a multi-use device which is capable of exceeding FDA Pre-enactment acoustic output (spatial peak-temporal average) intensity limits for fetal applications. The interaction of sound energy with tissue at sufficiently high levels and/or longer duration can produce biological effects (aka bioeffects) of either a mechanical or thermal nature.



DANGER

It is prudent to conduct an exam with the minimum amount and duration of acoustic output necessary to optimize the image's diagnostic value. The interaction of sound energy with tissue at sufficiently high levels can produce biological effects (aka bioeffects) of either a mechanical or thermal nature.

Concerns Surrounding Fetal Exposure

Always be aware of the acoustic output level by observing the Acoustic Output Display. In addition, become thoroughly familiar with the Acoustic Output Display and equipment controls affecting output.

Training

It is recommended that all users receive proper training in fetal Doppler applications before performing them in a clinical setting. Please contact a local sales representative for training assistance.

OB Type change

The ultrasound system includes measurements for the following studies: USA, Europe, Tokyo, Osaka and ASUM.

Select OB Type in **Utility > System > System Measure**.

Performing an Exam

NOTE

ASUM studies include the following measurements:

- ASUM: AC, BPD, and CRL
- ASUM 2001: AC, BPD, CRL, FL, HC, HL, and OFD

OB Report

The OB Report lists patient information, and all measurement and calculation data.

To view the OB Report, press the **Report** key on the control panel.

The OB Report has four sections of information:

1. Patient data
2. Measurement information
3. Calculation information
4. Report menu list

Patient Data

The Patient data section, at the top of the worksheet, lists information from the Patient Data Entry screen.

The following fields can be selected:

- FetusNo – If this is a multi-gestational patient, select the fetus in this field. **Fetus** selection can be adjusted to change the fetus.
- CUA/AUA – Select the ultrasound age calculation method
 - Composite Ultrasound Age (CUA) – Regression calculation
 - Average Ultrasound Age (AUA) – An arithmetic average

Select the method in this field or adjust the **Select CUA/AUA** control.

NOTE

CUA/AUA is only available when USA OB Type is selected in the Utility > System > System Measure menu.

Information can be entered in the following fields:

- FetusPos – Enter information about the fetus position.
- PLAC – Enter information about the placenta.

Measurement information

This section lists the results of all measurements.

- CUA or AUA - If this field is checked, the system uses the measurement to calculate the ultrasound age.
- Value - The measured value. If more than one measurement was made for an item, the system uses the specified method (average, maximum, minimum, or last) to determine this value.
- m1-m3 - Up to three measurement values for each item. If you make more than three measurements, the report uses the last three.

- Method - When there is more than one measurement for an item, this specifies the method used to calculate the measurement value listed in the Value column. Choices are average, maximum, minimum, last, or off. To change the method:
 1. Move the **Trackball** to the Method field.
 2. Press **Set**.
 3. Move the **Trackball** to select from the list.
 4. Press **Set**.
- AGE - The fetal age for this measurement.
- Range - The typical range of fetal age for this measurement.

Calculation Information

This section of the worksheet provides calculation choices and lists calculation results.

- EFW – Lists the parameters used to calculate EFW. This is followed by the calculation result.

To change which parameters are used:

- Select this field or press **Select EFW**.
- Select the desired parameters.
- EFW GP – Lists the source used to calculate EFW–GP (growth percentile). This is followed by the growth percentile.

To change the source:

- Select this field or press **Select GP**.
- Select the desired source.

The remaining calculation information shows ratios for several measurements, and the Cephalic Index (CI).

The worksheet shows if any of the ratios are out of range (OOR). Out of range indicates one of the following:

- The measurement is out of the normal range based on the gestational age that is calculated from the LMP. The system determines OOR from the ultrasound age compared to the gestational age. The gestational age is calculated from the last menstrual period or the estimated delivery date.
- The measurement is outside of the range for the data used in the calculation. That means that the measurement is either less than or more than the range of measurements used to determine fetal age based on the measurement.

Report Menu list

Use the report menu list to switch between different summary reports.

Performing an Exam

Anatomical Survey

The Anatomical Survey page provides a checklist that indicates which anatomy was imaged and its appearance.

Figure 2-15 Anatomical Survey

Origin LMP LMP BBT GA EDD(LMP)
Fetus=A/1 FetusPos PLAC

1 Normal 2 Clear 3 Save 4 Reset

ANATOMY	APPEARANCE	ANATOMY	APPEARANCE
HEAD		FOUR CHAMBER VIEW	
SPINE		STOMACH	
KIDNEYS		CORD INSERTION	
BLADDER		UPPER EXTREMITIES	
LOWER EXTREMITIES			

BIOPHYSICAL PROFILE	SCORE	COMMENTS
MOVEMENT		
TONE		
BREATHING		
FLUID		
REACTIVE NST		

Exit

1. **Normal/Seen** - Press to select normal/seen for all tables under APPEARANCE.
2. **Clear** - Clear all values in APPEARANCE tables.
3. **Save** - User can edit and save the current table as a default template.
4. **Reset** - Reset to factory default settings.

Editing

To activate the Anatomical Survey,

1. Select **Anatomy** on the OB Report screen.
2. Fill the required field.

Table 2-7 Anatomical Survey

Field	Description
FetusPos	Indicate the fetal position within the uterus.
PLAC	Identify the location of the placenta.

Field	Description
ANATOMY	1. Display the anatomy for each part. 2. User can add up to 11 additional anatomy. 3. User can edit each parameter in the anatomy column.
APPEARANCE	Indicate whether the appearance was seen, not seen, normal or abnormal.
BIOPHYSICAL PROFILE	The score is _ of 10 possible total points, depending upon the number of parameters entered. Enter the following information to assess the fetus's biophysical well-being.
MOVEMENT	Type/select 0, 1 or 2
TONE	Type/select 0, 1 or 2
BREATHING	Type/select 0, 1 or 2
FLUID	Type/select 0, 1 or 2
REACTIVE NST	Type/select 0, 1 or 2
COMMENTS	Free text

Select **Exit** to return to the Scan screen.

NOTE

The patient specific contents input on the Anatomical Survey page are returned to the factory default settings after starting a new patient.

OB Graphs

OB Graphs allow assessment of fetal growth compared to a normal growth curve. When a patient has completed two or more ultrasound exams, the graphs can also be used to look at fetal trending. For multi-gestational patients all fetuses can be plotted and the growth can be compared on the graphs.

The ultrasound system provides the following two basic types of graphs:

- **Fetal Growth Curve graphs** – show one measurement per graph. These graphs show the normal growth curve, positive and negative standard deviations or applicable percentiles, and ultrasound age of the fetus using the current measurement. For multi-gestational pregnancies, all fetuses can be viewed. If previous exam data is available, the graph can show fetal trending.
- **Fetal Growth Bar graph** – shows the ultrasound age and the gestational age based on patient data. Plots all measurements on one graph.

To View OB Graphs

To view OB graphs:

1. Press **Measure**.
2. Select **Graph**.

The system displays the OB Graph keys on the Touch Panel.

Performing an Exam

To Edit Patient Data

When working with graphs, the following patient data can be changed or entered:

- GA(LMP) – This field is computed using the LMP date on the Patient Data Entry screen. To change this field:

NOTE

This field can only be changed on the Fetal Growth Curve Graph in single view.

1. Move the **Trackball** to the field, which is left of the graph. To select the field, press **Set**.
The system displays a window with the GA weeks and days.
2. To select each field, move the **Trackball** to the field and press **Set**.
3. Type the correct weeks or days.
4. Select OK.

The system makes the following changes:

- GA (LMP) is now GA (GA) and shows the age you entered.
- In the Patient Data section, the GA changes.
- In the Patient Data section, The EDD (LMP) changes to EDD(GA) and shows an updated date, using the GA you entered.

The LMP is erased.

- FetusPos – Type information about the fetus position.
- PLAC – Type information about the placenta.

To Return from a Graph to the Scan Display

After viewing graphs, to return to the scan display, do one of the following:

- On the graph display, select Exit.
- On the Touch Panel, select **Graph**.

OB-Multigestational

Multiple Fetus

The ultrasound system allows measurement and reporting of the development of multiple fetuses. The system can report a maximum of four fetuses.

To Enter the Number of Fetuses

If more than one fetus is imaged during the exam, enter the number of fetuses in the Patient Data Entry Menu.

When an OB exam is started, the system automatically fills in the **Fetus #** field with 1. To change the number:

1. Move the cursor to the fetus number and press **Set** twice.

The number is highlighted.

2. Type the correct number and press **Set**.

The system displays a message to confirm changing the fetus number.

3. Select **Yes**.

To Identify Each Fetus

For measurements, calculations, and worksheet displays, the system labels each fetus A, B, C, or D. Each fetus is identified by a letter and the total number of fetuses. For example, fetus A/3 is fetus A from a total of 3.

When scanning, information can be entered about the fetus position and placenta location. The information can be entered in the Patient Data section of the worksheets and the graphs. Up to 23 characters can be entered in the **FetusPos** and **PLAC** fields.

Figure 2-16 Patient Data section of the OB Report

Origin	LMP	LMP	BBT	GA	EDD(LMP)
	Fetus=A/4	FetusPos	PLAC	Styles	Single

To Select a Fetus

During measurements and calculations, to change between fetuses, do one of the following:

- Adjust the **Fetus** selection.
- Move the **Trackball** to the Summary Window and select the fetus.

The fetus selection can be changed at any time during the exam.

NOTE

After changing to the next fetus, any measurements made are recorded and reported to that fetus. If any active measurement or calculation is not completed when changing the fetus, the system cancels the measurement or calculation.

To View Multiple Fetuses Data on Graphs

Gestational data for multiple fetuses can be displayed on fetal growth curve graphs. After making measurements for each fetus, select **Graph**.

1. To view the graph for each fetus, do one of the following:
 - Adjust the **Fetus** selection.
 - In the Patient Data section, move the **Trackball** to highlight the FetusNo field. In the list of fetuses, move the **Trackball** to select the desired fetus and press **Set**.
2. To display data for multiple fetuses on the same graph, select **Fetus Compare**.
The legend at the bottom of the graph shows the symbols and colors that represent each fetus.

To Show Fetal Trending for Multiple Fetuses

When there is data for more than one exam, fetal trending for multiple fetuses can be viewed and compared on one graph.

To view fetal trending for multiple fetuses:

1. Select **Graph**.
2. Select **Fetus Compare**.

3. Select **Plot Both**.

NOTE

Fetal trending for multiple fetuses can only be viewed in single graph display.

The symbol key for fetal trending and multiple fetuses is shown at the bottom of the graph.

To Compare Multiple Fetus Data on a Worksheet

Measurements of multiple fetuses can be listed and compared on the worksheet.

Select **Worksheet**, then select **Fetus Compare**.

When selecting **Fetus Compare**, the system lists the measurement results for each fetus on the Worksheet.

SonoBiometry (AFB)

SonoBiometry (AFB) is an alternative to the common fetal biometry measurements. It provides system suggested measurements for AC, BPD, FL, HL and HC which need to be confirmed by the user or can be changed manually.

Table 2-8 List of Abbreviations


Abbreviation	Full Form	Definition
AC	Abdominal Circumference	Abdominal Circumference of the baby's abdomen.
BPD	Biparietal Diameter	Transverse distance between the protuberances of the two parietal bones of the skull.
FL	Femur Length	The length in centimeters of the developing baby's femur (the femur is the long bone in the human arm).
HC	Head Circumference	Head Circumference of the baby's head.
HL	Humerus Length	The length in centimeters of the developing baby's humerus (the humerus is the long bone in the human arm).

Using SonoBiometry (AFB):

1. Press **Utility** and then select **System**.
2. Select **System Measure**.
3. Check the desired boxes under **SonoBiometry Option Selection** to enable system suggested measurements.
4. Select **OB**.
5. Press **Measure**.

6. Select BPD/HC/AC/FL/HL.

Figure 2-17 SonoBiometry (AFB) option



NOTE

SonoBiometry (AFB) measurements are supported on only recalled raw DICOM images and frozen images.

You can not use SonoBiometry (AFB) measurements on the following images:

- Images displayed in Dual/Quad view. In these views, an SonoBiometry (AFB) measurement will change to manual measurement automatically.
- Zoomed image.
- Rotated image.

NOTE

If a particular scanned image does not have valid fetus data, a warning message appears: "Auto measurement failure. Adjust calipers manually ...".

NOTE

If a particular image is not supported for SonoBiometry (AFB) measurement, a warning message appears: "Auto measurements do not operate on DICOM images. Using manual measurement ...".

GYN

The Gynecology exam category includes the following three studies:

- Generic. This study is common to all exam categories.
- General Gynecology. This study includes uterine, ovarian, ovarian follicle, and endometrium measurements.
- OB/GYN Vessel. This study includes the following vessels: uterine, ovarian, umbilical, middle cerebral artery, aorta, placenta, and descending aorta.

NOTE

The calculation formulas are listed in the Advanced Reference Manual.

B-Mode Measurements

In B-Mode, you make the measurements in the General Gynecology study. These measurements include:

- Uterine length, width, and height
- Ovarian length, width, and height
- Follicle
- Cervix
- Cervix Trace
- Endometrium

Follicle Measurements

Left and right ovary follicle measurements can be made from one, two, or three distances.

To select the left or right, rotate the rotary key on the control panel.

Endometrium Thickness

To measure the endometrium thickness, make one distance measurement.

1. Select **Endometrium**; an active caliper displays.
2. To position the active caliper at the start point, move the **Trackball**.
3. To fix the start point, press **Set**.
The system fixes the first caliper and displays a second active caliper.
4. To position the second active caliper at the end point, move the **Trackball**.
A dotted line connects the measurement points.
5. To complete the measurement, press **Set**.

The system displays the endometrium thickness in the Results Window.

Ovary Length, Width and Height

Length, width, and height of the left and right ovaries can be measured. Each measurement is a typical distance measurement made in the appropriate scan plane.

Typically, length and height are measured on the sagittal plane while the width is measured on the axial/transverse plane.

To measure ovarian length, width, or height:

1. Scan the patient's right or left ovary in the appropriate plane.
2. To select left or right, adjust the **Side** selection.
3. Select the **OV** folder, then select **OV L**, **OV W**, or **OV H**.
4. Perform a standard distance measurement.
 1. To position the active caliper at the start point, move the **Trackball**.
 2. To fix the start point, press **Set**.

The system fixes the first caliper and displays a second active caliper, if preset accordingly.

3. To position the second active caliper at the end point, move the **Trackball**.

A dotted line connects the measurement points.

4. To complete the measurement, press **Set**.

The system displays the distance value in the Results Window. After the first and second measurement, the system displays an active caliper for the next measurement.

5. To make the second and third distance measurements, repeat steps 3–4.

After you complete the length, width, and height measurements, the system displays the ovarian volume in the Results Window.

Uterus Length, Width and Height

Each of these is a standard distance measurement. Typically, length and height are measured on the sagittal plane while the width is measured on the axial/transverse plane.

To measure uterus length, width, or height:

1. Scan the patient in the appropriate scan plane.
2. Select the **UT** folder, then select **UT L**, **UT W**, or **UT H**.

An active caliper displays.

3. Perform a standard distance measurement.

1. To position the active caliper at the start point, move the **Trackball**.

2. To fix the start point, press **Set**.

The system fixes the first caliper and displays a second active caliper.

3. To position the second active caliper at the end point, move the **Trackball**.

A dotted line connects the measurement points.

4. To complete the measurement, press **Set**.

The system displays the distance value in the Results Window. After the first and second measurement, the system displays an active caliper for the next measurement.

4. To make the second and third distance measurement, repeat steps 2–3.

After you complete the third distance measurement, the system displays the uterine volume in the Results Window.

Cervix Measurements

Cervix measurements can be made from one distance or spline trace.

One Distance

1. Select **CX**; an active caliper displays.
2. To position the active caliper at the start point, move the **Trackball**.
3. To fix the start point, press **Set**.

The system fixes the first caliper and displays a second active caliper.

Performing an Exam

4. To position the second active caliper at the end point, move the **Trackball**.
A dotted line connects the measurement points.
5. To complete the measurement, press **Set**.
6. The system displays the cervix measurement in the Result Window.

Spline Trace

1. Select **CX Trace**; an active caliper displays.
2. To position the caliper at the start point, move the **Trackball**.
3. To fix the trace start point, press **Set**. The first caliper turns yellow. The second caliper appears at the same position as the first caliper and is green.
4. To position the second caliper, move the **Trackball** and press **Set**. The third caliper appears at the same position.
The spline trace requires at least three points to draw the trace. Continue setting the points of the trace until the desired points are set.
5. Press **Set** again after the last caliper is fixed to finalize the spline trace. All points are removed from the line and the spline trace turns yellow.

NOTE

Pressing **Set** twice finishes the trace measurement.

6. The system displays the cervix measurement in the Result Window.

Cardiac

Cardiology measurements offer two different types of measurement studies, Generic and Cardiac.

- Generic – Each exam category has a Generic study. The Generic studies provide you quick access to measurements.
- Cardiac – This study includes all cardiac measurements.

Naming Format for Cardiac Measurements

When making a measurement, select the abbreviation for the measurement on the Touch Panel. Most abbreviations are made using acronyms. The following table lists acronyms used for naming cardiac measurements.

Table 2-9 Cardiology Abbreviations

Acronym	Name
% STIVS	% Interventricular Shortening
A	Area
Acc	Acceleration
AccT	Flow Acceleration Time
ALS	Aortic Leaflet Separation
Ann	Annulus
Ao	Aorta

Acronym	Name
AR	Aortic Regurgitation
Asc	Ascending
ASD	Atrial Septal Defect
AV	Aortic Valve
AV Cusp	Aortic Valve Cusp Separation
AVA	Aortic Valve Area
AV-A	Aortic Valve Area by Continuity Equation
BSA	Body Surface Area
CI	Cardiac Index
CO	Cardiac Output
d	Diastolic
D	Diameter
Dec	Deceleration
DecT	Deceleration Time
Desc	Descending
Dur	Duration
EdV	End Diastolic Volume
EF	Ejection Fraction
EPSS	E-Point-to-Septum Separation
EsV	End Systolic Volume
ET	Ejection Time
FS	Fractional Shortening
FV	Flow Volume
FVI	Flow Velocity Integral
HR	Heart Rate
IVRT	IsoVolumetric Relaxation Time
IVS	Interventricular Septum
L	Length
LA	Left Atrium
LAA	Left Atrium Area
LAD	Left Atrium Diameter
LPA	Left Pulmonary Artery
LV	Left Ventricle

Performing an Exam

Acronym	Name
LVA	Left Ventricular Area
LVID	Left Ventricle Internal Diameter
LVL	Left Ventricle Length
LVM	Left Ventricular Mass
LVPW	Left Ventricle Posterior Wall
ML	Medial to Lateral
MPA	Main Pulmonary Artery
MR	Mitral Regurgitation
MV	Mitral Valve
MVcf	Mean Velocity Circumferential Fiber Shortening
MVO	Mitral Valve Orifice
OT	Outflow Tract
P	Papillary Muscles
PA	Pulmonary Artery
PAP	Pulmonary Artery Pressure
PDA	Patent Ductus Arteriosus
PEP	Pre-Ejection Period
PFO	Patent Foramen Ovale
PG	Pressure Gradient
PHT	Pressure Half Time
PI	Pulmonary Insufficiency
PISA	Proximal Isovelocity Surface Area
PR	Pulmonic Regurgitation
PV	Pulmonic Valve
PV-A	Pulmonic Valve Area by Continuity Equation
PVein	Pulmonary Vein
PW	Posterior Wall
Qp	Pulmonic Flow or CO
Qs	Systemic Flow or CO
RA	Right Atrium
RAA	Right Atrium Area
Rad	Radius
RAD	Right Atrium Diameter

Acronym	Name
RPA	Right Pulmonary Artery
RV	Right Ventricle
RVA	Right Ventricle Area
RVAW	Right Ventricle Anterior Wall
RVD	Right Ventricle Diameter
RVID	Right Ventricle Internal Diameter
RVL	Right Ventricle Length
RVOT	Right Ventricle Outflow Tract
s	Systolic
SI	Stroke Index
ST	Shortening
SV	Stroke Volume
SVI	Stroke Volume Index
T	Time
TA	Tricuspid Annulus
TAML	Tricuspid Annulus Medial to Lateral
TR	Tricuspid Regurgitation
TV	Tricuspid Valve
TVA	Tricuspid Valve Area
Vcf	Velocity Circumferential Fiber Shortening
Vel	Velocity
VET	Valve Ejection Time
Vmax	Maximum Velocity
Vmean	Mean Velocity
VSD	Ventricular Septal Defect
VTI	Velocity Time Integral

In this manual, the abbreviation for each measurement is listed in parenthesis after the measurement, as follows:

- Aortic Root Diameter (**Ao Diam**)
- Left Ventricle Posterior Wall Thickness, Diastolic (**LVPWd**)

For example, to measure the Aortic Root Diameter, select **Ao Diam** on the Touch Panel.

Auto EF Measurements

Automated Ejection Fraction (Auto EF) is a semi-automatic measurement tool used for measurement of the global EF (Ejection fraction). The Auto EF tool is used as an optional decision support tool.

The Auto EF tool tracks and calculates the myocardial tissue deformation based on feature tracking on B-Mode cine loops.

Auto EF is performed on either one or both apical 4-chamber or 2-chamber views, in any order.

The result is presented as Ejection Fraction value for each view and average Ejection Fraction for the whole LV. All values are stored to the worksheet after the results are approved.

NOTE

The Auto EF tool is intended for pediatric cardiology and is intended to be used in pediatric cardiology.

Acquisition

1. Create an exam, connect the ECG device and make sure to obtain a stable ECG trace.
2. Acquire B-Mode cineloops of an Apical 4 chamber view (4-ch) and an Apical 2 chamber view (2-ch).
 - The frame rate should be between 37 and 80 frames per second. A higher frame rate is recommended for high heart rate.
 - The Versana Premier/Versana Premier Lotus should be configured to store at least 100 ms before and after each heart cycle.
 - If the acquisition has more than one heart cycle, the analysis will be done on the second last heart cycle.
 - The entire myocardium should be visible.
 - A depth range that includes the entire left ventricle should be used

NOTE

The Auto EF processed image loop runs slower than the actual cardiac motion. To see the loop in correct playback speed, exit Auto EF.

Exit Auto EF

Press **Exit** to end Auto EF. Pressing Report, Patient or Scan also closes the Auto EF package.

NOTE

Do not disable any sub-measurement of **Auto EF** in **Utility > Measure**, as this will disable the whole set of Auto EF measurements.


NOTE

For consistent results, do not delete individual Auto EF measurements from the worksheet. Always delete a full column.

Vascular

Vascular measurements offer several different types of measurement studies:

Figure 2-18 Vascular Exam Category Touch Panel



- Generic – Common to all applications.
- Carotid
- LEA (Lower Extremity Artery)
- LEV (Lower Extremity Vein)
- TCD (Trans Cranial Doppler)
- UEA (Upper Extremity Artery)
- UEV (Upper Extremity Vein)
- Renal
- Bypass Graft
- UEV (Upper Extremity Vein) Map
- LEV (Lower Extremity Vein) Map

A vascular study is a group of particular vessels. Vessel exam calcs can be customized in the configuration menu.

When using Auto Vascular calculation, use the vessel keys on the Touch Panel to post assign vascular calculations. When not using Auto Vascular calculation, the vessel key is used for manual measurement.

Naming Format for Vessels

When you want to measure a vessel, on the Touch Panel you select the folder for the vessel. Many vessel folders are labeled with an abbreviation. The following table lists abbreviations used for naming vascular vessels.

Table 2-10 Vascular Vessel Abbreviations

Acronym	Name
ACA	Anterior Cerebral Artery

Performing an Exam

Acronym	Name
AComA	Anterior Communicating Artery
ATA	Anterior Tibial Artery
ATV	Anterior Tibial Vein
Axill	Axillary Artery
Axill V	Axillary Vein
BA	Basilar Artery or Brachial Artery
Ba V	Basilic Vein
Br V	Brachial Vein
CCA	Common Carotid Artery
Ceph V	Cephalic Vein
CFV	Common Femoral Vein
CHA	Common Hepatic Artery
Com Femoral	Common Femoral Artery
CIA	Common Iliac Artery
CIV	Common Iliac Vein
Com Iliac A	Common Iliac Artery
DFA	Deep Femoral Artery
DFV	Deep Femoral Vein
Dors Pedis	Dorsalis Pedis
DPA	Dorsalis Pedis Artery
ECA	External Carotid Artery
EIA	External Iliac Artery
EIV	External Iliac Vein
FV	Femoral Vein
GSV	Greater Saphenous Vein
ICA	Internal Carotid Artery (Transcranial Doppler)
ICA	Interior Carotid Artery (Carotid Artery)
IJV	Internal Jugular Vein
IMA	Inferior Mesenteric Artery
Inn	Innominate
IVC	Inferior Vena Cava
LSV	Lesser Saphenous Vein
MCA	Middle Cerebral Artery

Acronym	Name
Mcub V	Median Cubital Vein
Mid Hep V	Middle Hepatic Vein
MRA	Main Renal Artery
PCA	Posterior Cerebral Artery
PCoMA	Posterior Communicating Artery
Peron	Peroneal
POP	Popliteal
PTA	Posterior Tibial Artery
PTV	Posterior Tibial Vein
RA	Renal or Radial Artery
RV	Renal or Radial Vein
SMA	Superior Mesenteric Artery
SMV	Superior Mesenteric Vein
SUBC	Subclavian Artery
SUBC V	Subclavian Vein
SFA	Superficial Femoral Artery
TCD	Transcranial Doppler
TIPS	Transjugular Intrahepatic PortalSystemic Shunt
UA	Ulnar Artery
VERT	Vertebral Artery

IMT Measurement

The average of the intima media thickness can be measured for use as the index of arterial sclerosis.

IMT can be measured both on the posterior and the anterior walls of the vessel.

NOTE

Due to the physical properties of ultrasound imaging, the posterior IMT measurement is generally more accurate than the anterior IMT measurement.

IMT Measurement - Auto


Auto IMT automatically measures the thickness of the Intima Media on the far and near vessel walls. Near Wall IMT is the distance between the trailing edges of the adventitia and intima; the Far Wall IMT is the distance between the leading edges of the adventitia and intima.

Set up the parameters you want to record on the worksheet on the **Utility > Measure > M&A** page while are in the Carotid application. Select **CCA/ICA/BIF**, then click **AutoIMT Far/**

Performing an Exam


AutoIMT Near to check the box under **Parameter** (Average, Max, Min, Standard Deviation, Points, or Distance).

Figure 2-19 Configuring Auto IMT



In the Vascular Carotid application, the Auto IMT measurement is available.

Figure 2-20 Auto IMT Touch Panel



The following controls are available.

Table 2-11 Auto IMT Touch Panel Description

Parameter	Description
Worksheet	Select to view the Worksheet
IMT Far	Select to begin the Far Field IMT measurement.
IMT Near	Select to begin the Near Field IMT measurement.

Parameter	Description
AP	Anterior Posterior
Trans	Transverse
Length/Offset Rotary	Push to save Length/Offset as a preset. -40/+40 Length. At zero, you can freely adjust the length, but only vertically. Press key to save value as default. Offset distance, -20 (Left) / +20 (Right)
Overall / IMT Trace Fit / Intima	Adjusts (remeasures) the IMT automatically measured by the system.
Rt / Lt Side	Select Left / Right Side.
Cursor Select	Allows you to update cursor placement.

To measure the IMT:

1. In the Carotid application, press **Freeze**, press **Measure**.
2. Position the cursor, then select **IMT Far**.
3. Use the **Trackball** to set the length.

Or


Use the **Length / Offset** control on the Touch Panel to set the length and offset distance. The Offset key controls how far away from the vertical line the measurement starts. Length is the length of the tool itself. If set to zero, you can adjust it anywhere on the image.

4. Press **Set**.

Either adjust the trace prior to pressing the Print key, or press the Print key to store the image which also saves the measurement to the Worksheet.

To adjust the trace, use the **Overall IMT Trace Fit Intima** control on the Touch Panel. The Trace fit (up/down) adjusts the inter luminal line whereas the overall (rotate) adjusts both IMT lines.

Figure 2-21 Example of Auto IMT Far Measurement



5. Position the cursor, then select IMT Near.
6. Use the **Trackball** to set the length.
Or
Use the **Length / Offset** control on the Touch Panel to set the length and offset distance.
7. Press **Set**. "Store image to accept IMT measurement" displays in the message area. If the traces fit both layers of the wall, approve the measurement by pressing the **Print** key to store the image.

To adjust the trace prior to pressing the Print key, use the IMT Trace Fit control on the Touch Panel. The measurement is saved to the Worksheet.

NOTE

Since the IMT measurements are semi-automatic, the operator has to approve the detection by visual inspection before storing the results in the worksheet and report.

Auto Vascular Calculation

Auto Vascular Calculation enables the ultrasound system to detect and identify a cardiac cycle. It allows you to assign measurements and calculations during live timeline imaging, while the image is frozen, or in CINE. Peak values are detected for venous flow.

During cardiac cycle detection, the system identifies the cardiac cycle using calipers, vertical bars, and/or highlighting of timeline data. Use of identifiers is based on measurements and calculations selected by an operator for the current application. The system may place calipers at early systolic peak, peak systole, minimum diastole and end diastole. Vertical

bars may also be placed to indicate the beginning and end of the cardiac cycle. The peak and/or mean trace may be highlighted. You can edit the cardiac cycle identified by the system or select a different cardiac cycle.

Calculations to be displayed in the M&A Result window can be selected during live scanning or on a frozen image. These calculations are displayed at the top of M&A Result Window located adjacent to the image. These calculations are configurable by application, which means they can be set up the default calculations to be displayed for each application.

Activating Auto Vascular Calculation

To activate Auto Vascular Calculation, select **Auto Calcs** from **Live** (calculations displayed on the real-time image), or **Frozen** (calculations displayed on the frozen image).

To deactivate Auto Vascular calculation, select **Off**.

Set up Auto Vascular Calculation Parameters

- **Select Auto Trace**
To set continuous auto trace of the max or mean velocities:
 - Select Max or Mean using the **Trace Method** Touch Panel pull-down menu.
- **Select Trace Direction**
Trace Direction allows peak timeline data above, below, or composite (above and below) the baseline.
 - Select Positive, Negative or Both to set the peak timeline data.
- **Modify Calculation**
 1. Select the **Modify Calcs** Touch Panel key.
The Modify Calculation menu is displayed.
 2. Select which measurements and calculations are to be displayed in the Auto Vascular calculation window.

Select from the following parameters: PS, ED, MD, HR, TAMAX, PI, RI, Accel, PS/ED, ED/PS, AT, Volume Flow, PV.


Auto Vascular Calculation Exam

1. Preset the system.
2. Perform the scan and press **Freeze**.
3. Activate Auto Vascular Calculation.

Performing an Exam

The system performs a calculation automatically.

Figure 2-22 Auto Vascular Calculation




The Auto Vascular calculation is assigned to particular vessel measurements.

1. Press Measure to display the Measurement menu.
2. Select the location of the vessel (Prox, Mid, or Dist) and Side (Right or Left).
3. Select the desired vessel name from the Touch Panel.

Selected vessel measurements are automatically assigned with the Auto Vascular calculation. The results are then displayed in the Results Window.

Figure 2-23 Assigned Vessel



NOTE

To cancel the assignment, use the **Cancel Transfer** Touch Panel key.

During the course of an exam, the cardiac cycle may be indicated between two yellow bars; the peak trace and the mean trace may appear in green; calculation indicators appear on the spectral trace as a caliper identifier (these vary, depending on the selected calculation in the Results Window).

The right-most, most complete cycle is typically chosen to be the selected cardiac cycle. To select a different cardiac cycle:

- Move through CINE memory with the Trackball until the desired cardiac cycle is selected by the system.

NOTE

Several good cycles are required in front of the new cardiac cycle for this to be successful. Oftentimes, this is problematic near a freeze bar.

- Use the **Cycle Select** control to cycle to a different cardiac cycle.

To move the systole or diastole position:

- Use the **Cursor Select** control to move the start systole position or the end diastole position.

Manual Vascular Calculation

Calculations can be performed manually when Auto Doppler Calculation is not activated.

Performing an Exam


1. Press **Measure**.

If necessary, you can select another Exam Calc and then select parameters from Modify Calculation.

2. Select the location of the vessel (Prox, Mid, or Dist) and Side (Right or Left).
3. Select the desired vessel folder.

The Measurement menu is displayed.

Figure 2-24 Measurement Menu Example



4. Make the required measurements according to the system, or select your preferred measurements.

To select vascular measurements

Your system is set up to show the measurements that you usually make for each vessel. To make a measurement that is not shown for the selected vessel:

1. Enter **Utility > Measure**, select **Vascular**.
2. The system displays all possible vessel measurements.
3. Select **Doppler** tab.
4. Select the desired measurement.

NOTE

The following instructions assume that you first scan the patient and then press Freeze.

Urology

Urology measurements offer three different types of measurement studies:

- Generic—Common to all applications.
- Urology

- This chapter describes Urology B-Mode measurements.
- The Urology M-Mode measurements are common to other applications.
- The Urology Doppler measurements are common to other applications.
- Pelvic Floor.

NOTE

Bladder(0.7) Vol, Bladder Vol, Post Void Vol, Prostate Vol, Renal Vol, Renal (0.8) Vol and Volume can be displayed on the Touch Panel if preset at the Utility > Measure screen.

Urology B-Mode Measurements

In B-Mode, the Generic Exam Calcs for Urology includes the following measurements:

- % Stenosis
- Volume
- Angle
- A/B Ratio
- Elasto
- E Ratio

The following measurements are located specifically in the Urology Exam Calcs. Those specific measurements (Bladder Volume, Prostate Volume and Renal Volume) are listed on the following pages.

Press **Measure** key and select **Exam Category** on Touch Panel. Then the following exam category is displayed.

NOTE

Bladder(0.7), Bladder(0.5), Post Void, Prostate, Renal, Renal (0.49), Volume and STVOL can be displayed if preset at the **Utility > Measure** screen.

Bladder Volume

This calculation uses a standard distance measurement. Length is typically measured in the sagittal plane. Width and height are measured in the axial plane.

To measure Bladder Volume:

1. Scan the patient in the appropriate scan plane.
2. Select the **Bladder** folder, an active caliper displays.
3. Perform a standard distance measurement.
The system displays the distance value in the Results Window.
4. To make the second and third distance measurement, repeat steps 2–3.

After you complete the third distance measurement, the system displays the bladder volume in the Results Window.

Renal Volume

This calculation uses a standard distance measurement. Length is typically measured in the sagittal plane. Width and height are measured in the axial plane.

To measure Renal Volume:

Performing an Exam

1. Scan the patient in the appropriate scan plane.
2. Select **Rt** or **Lt** side.
3. Select the **Renal** folder, an active caliper displays.
4. Perform a standard distance measurement.
The system displays the distance value in the Results Window.
5. To make the second and third distance measurement, repeat steps 2–4.

After you complete the third distance measurement, the system displays the renal volume in the Results Window.

Prostate Volume

This calculation uses a standard distance measurement. Length is typically measured in the sagittal plane. Width and height are measured in the axial plane.

To measure Prostate Volume:

1. Scan the patient in the appropriate scan plane.
2. Select the **Prostate** folder, an active caliper displays.
3. Perform a standard distance measurement.
The system displays the distance value in the Results Window.
4. To make the second and third distance measurement, repeat steps 2–3.

After you complete the third distance measurement, the system displays the prostate volume in the Results Window.

PSA Measurement

Entering the value of PSA (Prostatic Specific Antigen) and PPSA Coefficient at the Urology Patient screen, automatically calculates PSAD and PPSA.

The values are displayed on the Worksheet and Report (if set appropriately on the Report Designer page).

Figure 2-25 Urology Patient Screen

ABD	OB	GYN	CARD	VAS	UR	SMP	PED	MSK	THOR
PSA						0.000		ng/ml	
PPSA Coefficient 1						0.12			
PPSA Coefficient 2						0.12			
Operator						ADM			
Exam Description									
Accession #									
Perf.Physician									
Ref.Physician									
						Clear			

Figure 2-26 Measurement result window

1	Prostate L	6.12 cm
2	Prostate H	5.23 cm
3	Prostate W	7.72 cm
	Prostate Vol	129.21 ml
	PSAD	0.00
	PPSA(1)	15.51
	PPSA(2)	15.51
+	d	15.04 cm
	L	0.00 cm

PSAD: Prostatic Specific Antigen (PSA) Density – defined as: $PSAD = PSA/Volume$

PPSA: Predicted Prostate Specific Antigen – defined as: $PPSA = Volume \times PPSA \text{ Coefficient}$

Report


- For the prostate volume calculation, you can select the method "m1, m2..." in addition to Avg., Max., Min. and Last.
- The value of PSA and PPSA are displayed.

Performing an Exam

Pelvic Floor Measurements

Pelvic floor measurements can be performed in the Pelvic Floor study. The measurements are located in the Exam Calc folder in the Urology preset.

Figure 2-27 Pelvic Floor Study Touch Panel



BN (Bladder Neck) Rest

Obtain an image with the patient at rest (relaxed).

1. Create a straight line (zero or baseline) to line up with the inferior/posterior of symphysis pubis bone.
2. Once the baseline is positioned, a caliper appears. Position the caliper at the anterior margin of the bladder neck. A positive number displays since the caliper is placed below the baseline.
3. A distance is calculated in millimeters.

BN (Bladder Neck) Stress

Obtain an image after the patient performs the Valsalva maneuver.

1. Create a straight line (zero or baseline) to line up with the inferior/posterior of symphysis pubis bone.
2. Once the baseline is positioned, a caliper appears. Position the caliper at the anterior margin of the bladder neck.

If the bladder neck is below the baseline, the Bladder Neck Stress is a positive number. If the bladder neck is above the baseline (closer to the transducer face), the number is negative.

BN (Bladder Neck) Descent

The Bladder Neck Descent is a calculation that should be calculated after measuring the Bladder Neck Rest and Bladder Neck Stress.

BND = Bladder Neck Rest - Bladder Neck Stress

NOTE

If the Bladder Neck Stress is a negative number, it becomes positive and is added to the bladder neck rest measurement.

DWT (Detrusor Wall Thickening)

Three distance measurements of the bladder wall dome are calculated into a mean dimension and displayed in millimeters.

UT (Uterine) Descent Max

1. Create a straight line (zero or baseline) to line up with the inferior/posterior margin of symphysis pubis bone.
2. Measure using a 2-caliper dimension to the inferior position of the uterus in a stress image and display in millimeters

Rect Amp Des Max (Rectal Ampulla Descent Max)

1. Create a straight line (zero or baseline) to line up with the inferior/posterior margin of symphysis pubis bone.
2. Measure using a 2-caliper dimension to the inferior position of the rectal ampulla in a stress image and displayed in millimeters

Rectocele (Depth and Width)

Two 2-caliper diameter measurements to measure depth and width of the rectocele. Displayed in millimeters.

Lev Hiat Stress (Levator Hiatus Stress)

Two 2-caliper diameter measurements and calculate an area displayed as cm squared.

Residual Urine

Two 2-caliper diameter measurements calculate as:

(x) times (y) times 5.9 minus 14.9 equals Residual Volume displayed in ml.

Auto Bladder Volume Measurement (Whizz On)

Auto Bladder Volume automatically measures the Bladder Volume. It's main application is to measure three longest orthogonal lines (L, H and W) from two bladder slices.

It provides "Auto Bladder LxH" and "Auto Bladder W" to take L, H and W values, then automatically calculate the Bladder Volume.

Set up the parameters you want to record on the report in **Utility>Measure>M&A** page. Select the **Auto Bladder Volume Tool**.

The following controls are available.

Table 2-12 Auto Bladder Volume description

Parameter	Description
L	L value from bladder slice.

Parameter	Description
H	H value from bladder slice.
W	W value from the bladder slice.
Vol	Bladder volume calculated by $0.7 \times L \times H \times W$ or $0.5 \times L \times H \times W$.

To measure the bladder volume,

1. Select **Bladder** application.
2. Enter **Utility>Application>Whizz** page, check **Preset** as **Bladder** and configure **Auto IQ Optimization**, **Whizz CF**, **Auto Measurement** and **Auto Annotation**.
3. Press **Whizz** key on the control panel to turn on Whizz.
4. Activate the M&A package, press **Measure**.
5. **Auto Bladder(0.7)** is selected automatically. If user wants to select **Auto Bladder (0.5)**, press **Return** to continue.
6. Place the caliper anywhere outside the bladder and position it. Press **Set**.
7. The first caliper is anchored in place and rectangle appears, free to adjust.
Move the second caliper to the right or to the left to generate a dotted rectangle connecting the two cursors.
8. Press **Set** and complete the measurement. Annotation of Bladder will appear automatically in the screen.
9. After a measurement, the length of lines are displayed in the result window as "Bladder L/H/W and Bladder Vol".
10. User could select and adjust the calipers of lines by "Cursor Select".

NOTE

Do not select the region of interest including many non-bladder area. Do not select the region of interest inside the bladder area.

NOTE

The accuracy for auto bladder volume measurement is $\pm 30\%$.


Pediatrics

Pediatrics measurements offer two different types of measurement studies:

- Generic. The Generic Calculations study is common to all applications.
- Pediatric Hip (PedHip).
 - This chapter describes Pediatrics B-Mode measurements.
 - The Pediatrics M-Mode measurements and Doppler measurements are common to other applications.
 - The Pediatrics Doppler measurements are common to other applications.

Pediatrics Hip

Figure 2-28 Pediatrics B-Mode Measurement



Hip Dysplasia Measurement

The HIP calculation assists in assessing the development of the infant hip. In this calculation, three straight lines are superimposed on the image and aligned with the anatomical features. The two angles are computed, displayed, and can be used by the physician in making a diagnosis.


The three lines¹ are:

- 1 - The baseline connects the osseous acetabulum convexity to the point where the joint capsule and the perichondrium unite with the iliac bone.
- 2 - The inclination line connects the osseous convexity to labrum acetabulare.
- 3 - The Acetabulum roof line connects the lower edge of the os ilium to the osseous convexity.

The α (Alpha) angle is the supplement of the angle between 1 and 3. It characterizes the osseous convexity. The β (Beta) angle is the angle between lines 1 and 2. It characterizes the bone supplementing additional roofing by the cartilaginous convexity.

Figure 2-29 Hip Dysplasia

Performing an Exam



Anatomical Landmarks

- a** Ilium
- b** Iliac Bone
- c** Labrum
- d** Bony Roof
- e** Cartilaginous acetabular roof
- f** Femoral Head

¹Source: R GRAF, Journal of Pediatric Orthopedics, 4: 735-740(1984)

To make a Hip Dysplasia measurement:

1. From the Touch Panel, select either the **right** or **left side** (orientation) and then select **Beta Alpha HIP**.

A horizontal dotted line displays.

2. To place the baseline, move the **Trackball**. Position the crosshairs edge at the osseous convexity of the ilium.
3. To rotate or change inclination, adjust the **Ellipse** control or **Hip Rotate**.
4. To fix the baseline, press **Set**.

The system displays a second dotted line at an angle.

5. To place the line along the inclination line of the osseous convexity to labrum acetabulare, move the **Trackball**.
6. To rotate or change inclination, adjust the **Ellipse** control or **Hip Rotate**.
7. To fix the second measurement line, press **Set**.

The system displays a third dotted line at an angle.

8. To place the caliper along the acetabular roof line, move the **Trackball**.
9. To rotate or change inclination, adjust the **Ellipse** control or **Hip Rotate**.
10. To fix the third measurement line and complete measurement, press **Set**.

The system displays the hip measurements (α and β) in the Results Window.

Alpha HIP

The Alpha HIP measurement measures the angle between the iliac baseline and the bony roof line. To make an Alpha HIP measurement:

1. From the Touch Panel, select either the **right** or **left side** (orientation) and then select **Alpha HIP**.

A horizontal dotted line displays.

2. To place the baseline, move the **Trackball**. Position the crosshairs edge at the osseous convexity of the ilium.
3. To rotate or change inclination, adjust the **Ellipse** control or **Hip Rotate**.
4. To fix the baseline, press **Set**.

The system displays a second dotted line at an angle.

5. To place the caliper along the acetabular roof line, move the **Trackball**.
6. To rotate or change inclination, adjust the **Ellipse** control or **Hip Rotate**.
7. To fix the second measurement line, press **Set**.

The system displays the alpha hip measurement (α) in the Results Window.

d:D Ratio Measurement

The d:D Ratio measurement measures the percentage of the femoral head coverage under the bony roof. To make this measurement:

1. From the Touch Panel, select either the **right** or **left side** (orientation) and then select **d:D Ratio**.
A horizontal dotted line displays.
2. Use the **Trackball** to place the baseline along the ilium. Position the crosshairs edge at the osseous convexity of the ilium.
3. Use the **Ellipse** control to adjust or change inclination or **Hip Rotate**.
4. Press **Set** to fix the baseline.
5. The system displays a circle representing the femoral head. Use the **Trackball** to position the circle.
6. Use the **Ellipse** control to size the femoral head circumference.
7. Press **Set** to fix the femoral head circumference.

The system displays the d:D ratio for the femoral head in the Results Window.

Auto vs. Manual Calculations

The same calculations can be performed using either manual or auto calcs.

Manual Calcs

To perform manual calcs:

1. To turn Auto Calcs off and perform manual measurements, choose **Auto Calcs > OFF** on the PW tab of the Touch Panel.
2. After obtaining a waveform, press **Measure**. Choose the appropriate vessel folder or calculation. The system walks you through the measurement.

NOTE

To program which calculations are done under manual calcs when using measurement folders for measuring specific vessels, press the Utility key. Select Measure > Doppler and program your manual calcs (Auto Calcs OFF). Each vessel must be programmed individually and saved after each change.

Performing an Exam

Auto Calcs

To perform auto calcs:

1. Ensure that the auto calcs function is on by choosing **Auto Calcs > Frozen** or **Live** on the Doppler tab of the Touch Panel.
 - Live: Auto calculation activates when the system is in real-time.
 - Frozen: Auto calculation activates when you press Freeze.
 - Off
2. After obtaining a waveform, press **Measure**. Choose the appropriate vessel folder, side and location. The measurements that are pre-programmed are performed automatically and entered in the worksheet.

To modify auto calcs:

1. Select **Modify Auto Calcs** on the Touch Panel.
2. Choose the measurements to be performed with this preset.
3. To save these measurements:
 - If this is a temporary change, press Return.
 - If this is a permanent change, select **Save as default**.

The measurements are saved and can be performed with the auto calcs function.

Edit Auto Calcs

Auto Calcs can be edited after taking an Auto Trace measurement.

1. After taking an Auto Calc with a trace, select the measurement result on the result window. The Edit Trace menu window appears.

NOTE

If the system cannot take the trace data correctly from the image, Edit Trace does not work.

2. Select Edit Trace. The first caliper (manual trace caliper) appears on the center of the image. Use the **Trackball** to move the caliper on the trace line to the start point.

NOTE

To cancel Edit Trace at this time, press **Clear**, **Scan**, or **Freeze**.

3. Press **Set** to fix the first caliper. The second caliper appears. Edit the trace manually using the second caliper.

The Ellipse control is used to edit the trace.

NOTE

When pressing the **Clear** key once at this time, the second caliper disappears and the first caliper appears in the center of the image.

NOTE

If you press **Scan** or **Freeze** at this time, the caliper is automatically fixed and the result window updates.

4. Press **Set** to fix the second caliper. The trace and the result window are updated. The data is retaken from the trace and updated.

NOTE

While in Edit Trace, Cursor Select is disabled.

The trace data (TAMAX and TAMEAN) is updated, but the other selections (e.g. PS, ED) are not updated by trace. The points can be edited using **Cursor Select** if needed.

5. Repeat Edit Trace as needed.

Modify Auto Calcs


When you select this key, the Modify Calculation menu is displayed as below. In this menu, you select parameters to display in the Auto Vascular Calculation window. Only parameters that can be used by the calculation are displayed.

Select **Save as Default** to save the selected parameters as the default calculations for this application.

Select **Return** to return to the previous Touch Panel screen.

If you select **PV**, all selected parameters are turned off. When you deselect **PV**, the system returns to the previously selected calculation.

Figure 2-30 Modify Auto Calculation Menu



Wide Dual Screen Measurements

To perform measurements in wide dual screen, select measurement item from Touch Panel.

NOTE

Ensure that the "Dual Screen" has been selected in the Automatic Wide Screen... portion of the **Utility>System>System Imaging**.

Performing an Exam

1. Press **Dual Live** button on the touch panel to enter wide dual screen.
2. Press **Measure** on the control panel, the Result window displays at the left bottom of the screen. To perform measurement, select measurement item from left measurement list or select measurement item from Touch Panel. To switch from different mode measurement package, press Left/right arrow keys on the Alphanumeric keyboard.

NOTE

The location of Result Window and measurement item in wide dual screen is different from single screen and dual screen.

Viewing and Editing Worksheets

As you complete measurements, the system puts measurement data in the appropriate worksheets

NOTE

Worksheets are not saved if the system crashes.

For the same measurement item, the most recent 6 values can be averaged.

To view a worksheet

To view a worksheet, select **Report** on the Touch Panel.

The system displays the worksheet for the current study.

If a worksheet has more data on a second page, to view the next page, rotate the **Page Change** knob or select directly on the Touch Panel.

To return to scanning, do one of the following:

- Select **Report**.
- Press **Esc**.
- Select the **Exit** button by **Cursor**.

To edit a worksheet



HINTS

Some fields on the report are view only, and others you can change or select. To easily see which fields you can change or select, move the Trackball. As the cursor moves over a field that you can change or select, the field is highlighted.

Change Data

1. Move the **Trackball** to the field you want to change.
2. Press **Set**.
3. Type the new data in the field. The new data is displayed in blue to indicate that it was manually entered.

Delete or Exclude Data

1. Move the Trackball to the field you want to delete or exclude, press the Cursor. The field is highlighted.
2. Do one of the following:
 - To delete the field, select **Delete Value**.
 - To exclude the field, select **Exclude Value**.
The data in the field is not visible and is not included in worksheet calculations.
 - To include a value that you previously excluded, select **Exclude Value**.

Type a comment

1. Select **Comments**. The Comments window opens.
2. Type comments about the exam.
3. To close the Comments window, select other Menu.

Turn the volume measurement value off


To turn the volume measurement value off, select the method type **Off**.

Delete All Worksheet Values

You can delete all measurement values on a worksheet.

1. When the Worksheet is displayed on the monitor, press the **Clear** key, the following warning message appears:

Figure 2-31 Warning Message



2. Select **OK** to delete all.
Select **Cancel** to cancel the deletion.

Clinical Measurement Accuracy

Basic Measurements

The following information is intended to provide guidance to the user in determining the amount of variation or measurement error that should be considered when performing clinical measurements with this equipment. Error can be contributed by equipment limitations and improper user technique. Be sure to follow all measurement instructions and develop uniform

measurement techniques among all users to minimize the potential operator error. Also, in order to detect possible equipment malfunctions that could affect measurement accuracy, a quality assurance (QA) plan should be established for the equipment that includes routine accuracy checks with tissue mimicking phantoms.

Please be advised that all distance and Doppler related measurements through tissue are dependent upon the propagation velocity of sound within the tissue. The propagation velocity usually varies with the type of tissue, but an average velocity for soft tissue is assumed. This equipment is designed for, and the accuracy statements listed on are based on, an assumed average velocity of 1540 m/s. The percent accuracy when stated applies to the measurement obtained (not the full scale range). Where the accuracy is stated as a percent with a fixed value, the expected inaccuracy is the greater of the two.

Table 2-13 System Measurements and Accuracies

Measurement	Units	Useful Range	Accuracy	Limitations or Conditions
Depth	mm	Full Screen	≤10%	
Angle	degree	Full Screen	≤5%	
Distance:				
Axial	mm	Full Screen	<5%	
Lateral	mm	Full Screen	<5%	Linear Probes
Lateral	mm	Full Screen	<5%	Convex Probes
Lateral	mm	Full Screen	<5%	Sector Probes
Circumference:				
Trace	mm	Full Screen	≤10%	Linear Probes, Convex Probes, Sector Probes
Ellipse	mm	Full Screen	≤5%	Linear Probes, Convex Probes, Sector Probes
Area:				
Trace	mm ²	Full Screen	≤5%	Linear Probes, Convex Probes, Sector Probes
Ellipse	mm ²	Full Screen	≤5%	Linear Probes, Convex Probes, Sector Probes
Time	s	Timeline Display	<5%	M mode, AM mode, CM mode, PWD mode, CWD mode
Slope	mm/s	Timeline Display	≤10%	M mode, AM mode, CM mode
Doppler SV Position	mm	Full Screen	≤2 mm	PWD mode
Doppler Velocity	cm/s	Form 0 to 100 cm/s	<15%	PWD mode, CWD mode
		From 100 to 130 cm/s	<10%	

Measurement	Units	Useful Range	Accuracy	Limitations or Conditions
Doppler Angle Correction	cm/s	From 0-80°	≤5%	PWD mode

Chapter 3

After The Exam Is Over

In this section

Preset.	152
Electronic Documentation (Not applicable for China).	193
System Care and Maintenance.	196
Quality Assurance.	233
Image Quality Check (IQC).	248
Privacy and Security.	249
Contact Information.	250

Preset

Preset Menus provides the following functionality:

- **System.** View and update general system configuration settings, measurement settings, backup and restore data, Peripherals, User Configurable Key, User Defined Gesture, Probe Buttons and configuration files.
- **Imaging.** View and update exam and imaging parameters.
- **Comments.** Set up comment libraries by application.
- **Body Patterns.** Set up body pattern libraries by application.
- **Application.** Configure application- and user-specific settings.
- **Test Patterns.** Helps configure system settings.
- **Connectivity.** Define connection and communication setup, including exam dataflow information.
- **Measure.** Customize exam studies, create measurements, set up manual sequencing, and create OB Tables.
- **Whizz Report.** Allows you to manage the report template, import the report template, and export Whizz Report Editor for creating user customized report template.
- **Admin.** Perform system administrator activities such as setting up user IDs and logon formats.
- **Service.** Activates the Service Browser.
- **Scan Coach.** Create, import/export, and manage Scan Coach programs.
- **Scan Assistant.** Create, import/export, and manage Scan Assistant programs.
- **Imaging Preset Manager.** Manage Imaging Presets
- **Whizz.** Customize Whizz setting.
- **Search.** User can search for a parameter on the Utility pages (Measure, Reports, and Service pages cannot be searched.)
- **3D/4D.** Real-time 4D and Static 3D scanning.
- **LDAP**(Lightweight Directory Access Protocol). Allows you to enable/disable LDAP authentication, and configure LDAP.
- **Barcode.** Manage barcode configuration.

To access these functions, select the **Utility** on the Touch Panel, then select the appropriate menu key.

System Presets

System presets allows you to view or change the following parameters:

- **General** - Location, Date/Time, System Prompt Sound, Bi-plane Simultaneous Layout, Report File Name Composition, Title Bar, Trackball, Key Usage, Utility, Scan Assistant, Audio and Breast Care configuration
- **System Imaging** - Biopsy Guides, Follow Up Tool, Image Label Color, Image Label Layout, Auto Zoom Linear Probe Images at Shallow Depth..., Controls, Display, Automatic Wide Screen... and QAnalysis Statistic Value configuration
- **System Measure** - Measurement, Cursor and Results Window configuration, SonoBiometry Option Selection

- **Backup/Restore** - Backup, Media, EZMove, EZBackup, Restore, Detailed Restore of User Defined
- **Peripherals** - Print and Store Options, Removable Media, Standard Printer Properties, Default Printer and Setup configuration
- **User Configurable Key** - User Defined Key, Keyboard Key and User Defined Trackball Key configuration
- **User Defined Gesture** - User Defined Gesture configuration.
- **Probe buttons** - L4-20t probe buttons configuration.
- **About** - Software, Base Image, Patents, Scan Coach Addendum, System, System Image Patches and Additional Information.

Changing system parameters

To change system parameters:

1. Select **Utility** on the Touch Panel.
2. Then select **System**. The System screen displayed.
3. On the monitor display, move the **Trackball** to select the tab that has the information you want to change.
4. Select values for the parameters you want to change.
5. To save the changes, press **Save** button. Press **Exit** to return to scanning. In some cases, you may need to reboot the system for the change to take effect.

Keyboard Language Setup

1. Press **Utility > System > General**, select desired language under **Language (restart needed)** and **Input Language** drop-down menu.
2. Press **Save**.
3. System restart pop-up appears, press **Restart now** to restart the system.

NOTE

To have the settings take effect, you MUST restart the system.

NOTE

Press Alt+Shift to change the input language.

Imaging Presets


Imaging screens allow you to specify parameters for the following:

- B-Mode (B)
- Color Flow Mode (CF)
- Power Doppler Imaging (PDI)
- Elastography (ELASTO)
- M-Mode (M)

After The Exam Is Over

- Anatomical M-Mode (AMM)
- Pulse Wave Mode (PW)
- Continuous Wave Mode (CW)
- Harmonics (HAR)
- B Flow Mode (BF)
- B Flow Color (BFC)
- Ref
- Contrast (CON)
- Tissue Velocity Imaging (TVI)
- Tissue Velocity Doppler (TVD)
- General

Figure 3-1 Imaging Preset Example



1. Preset/application dependent setup parameters
2. Probe dependent setup parameters

Changing imaging presets

To change imaging presets:

1. On the Touch Panel, select **Utility**.
2. On the Touch Panel, select **Imaging**.
The system displays the Imaging screens.
3. In the row across the top of the screen, select the mode.

The system displays two sets of parameters and settings. The left column lists all settings for the exam (for example, Abdomen). The right column(s) list settings that apply only to the exam and probe combination.
4. In the Preset list, select the exam.
5. In the Probe list, select the probe.
6. To change a parameter, do one of the following:
 - Select the value from a list
 - Select one value from a choice of two or more buttons
 - Select or clear a check box
7. After changing the parameters, to save the changes, select the Save button.

NOTE

When you Save changes to imaging parameters, the system saves changes to all modes, not just the mode currently displayed.

If you have problems with imaging, you can return parameters back to the original settings. Select the exam, probe, and mode, and then select Reload Factory Defaults. The system returns the selected parameters to the original settings.

For information about the specific parameters, refer to Chapter 2 Optimizing the Image.

General

You can specify a default probe per application and a default application per probe, ECG Display or Sync Mode.

Default probe per application

1. To specify a default probe per application, select **Utility > Imaging > General**.
2. Check the parameter if you want to start it automatically.
3. Select the default probe from the pull-down menu.

Default mode and application per probe

1. To specify a default application per probe, select **Utility > Imaging > General**.
2. Under Probe, specify the desired mode and application from the pull-down menu.

Checkmark the following fields when you want the system to activate a certain display. Values vary by probe.

- Simultaneous
- ECG
 - ECG Display
 - Sync Mode

If your system has the following options, you must assign the option to the button per application and probe.

- B or HAR(Harmonics)
- Acoustic Output (%)

NOTE

If you use the user-defined preset as default preset, re-select the preset in Utility whenever you overwrite the user-defined preset.

Imaging Preset Manager

Selecting an Application Preset

The exam category preset that best describes the desired exam to be performed is chosen after the exam category is selected. The factory default preset selections are displayed on the menu area.

Use these parameters as a starting point for the exam.


Creating a User-Defined Application Preset

To create a User-Defined Application Preset,

After The Exam Is Over

1. Select the **Probe** tab on the Touch Panel.
2. Select the Model and Application you want to use as a basis for the new Application Preset.
3. Select **Create New**.

Figure 3-2 Application Preset Menu



4. The **Create New Application** menu appears. Enter the name and check the "**Add to My Preset**" box, then press **Create**.

NOTE


The new user application is based upon the current exam category and application, plus any modifications you have made, including the comment library and M&A calcs.

NOTE

The name of the new application cannot include spaces or symbols. However, the name can include numbers and letters.

- To view/edit the parameters for the user-defined preset, press **Utility > Imaging Preset Manager**.

Figure 3-3 Preset Manager



If you change the settings for this application, make sure to save the changes via **Save > Overwrite (user application)**.

NOTE

If you select Reset for the user-defined application that you created, the settings for this user-defined application revert back to the factory settings for the exam category and application it was based upon.

NOTE

If you use the user-defined preset as default preset, re-select the preset in Utility whenever you overwrite the user-defined preset.

NOTE

You can delete, rename, edit a user-defined application. And you can also create a new application preset by selecting **Create New**.

Arranging Presets Position

On this screen, you can specify which and where you want the application presets to display on the exam's Application Preset Menu.

To map the preset to the preset keys:


Select the mapping grid under **My Preset/Preset/Fifth Column** and select the preset you want to display under **Imaging Preset Manager**, then press >>. The selected preset will be under **My Preset/Preset/Fifth Column**.

OR

After The Exam Is Over

Select the preset under **Imaging Presets Manager** and drag to any position you want to locate under **My Preset/Preset/Fifth Column**.

Figure 3-4 Mapping the preset keys



User can move the location of the preset displaying on the Application Preset Menu via **Utility>Imaging Preset Manager**.

To reposition a preset,

1. Select the Preset you want to move under **My Preset/Preset/Fifth Column**.
2. Use the **Up/Down** button on the left side to move the preset to the next grid, or press the preset selected and drag to a new grid.

Fifth Column

Fifth Column locates on the right side of the Touch Panel screen. It contains **Frequent Presets**, **Recent Presets** and **Customized Presets**.

Frequent Presets - Most frequent used presets.

Recent Presets - Most recent used presets.

Customized Presets- User can use >> key to move or directly drag the presets to this page.


The Fifth Column on the Touch Panel provides quick access to some presets from touch panel. Swipe the fifth column with one finger from right to left/left to right to change between Frequent Presets, Recent Presets and Customized Presets.

NOTE

Swipe black space can help to prevent accidentally Hit/Select any button on the touch panel.

User can move the position of customized presets displaying on the Touch Panel via **Utility>Imaging Preset Manager** (on touch panel).

Figure 3-5 Fifth Column Configuration



To reposition a Preset on the Touch Panel grid,

1. Select the Preset you want to move under **Customized Presets**.
2. Use the **Up** or **Down** button to move the Preset to the next grid, or press and drag the selected Preset to any grid preferred.

Updating User Presets

You can edit, reset to factory default, or delete any user preset you create, as long as you have selected it in the "Available Imaging Presets" column on the left.

Editing imaging Parameters

To view/edit the parameters for the user-defined preset,

After The Exam Is Over

1. Adjust the image while in the user preset you want to edit.
2. Press **Probe** on the Touch Panel.
3. Select [**Preset Name**] and press **Overwrite**.
4. From the **Utility > Imaging Preset Manager**.

To view/edit the parameters for the user defined preset

1. Press **Utility > Imaging Preset Manager**. Select the user preset you want to view/edit from Available Imaging Presets Column.
2. Press **Edit**. The Imaging page appears.
3. Edit the presets as necessary and press **Save**.

OR

1. Adjust the image while in the user preset you want to edit.
2. Press **Probe** on Touch Panel.
3. Select [**Preset Name**] and press **Overwrite**.

Renaming a User Preset

To rename a user preset,

1. Press **Utility > Imaging Preset Manager**. Select the user preset you want to rename.
2. Press **Rename**. The Rename Preset pop-up menu appears.
3. Type the new name and press **Rename**.

Deleting a User Preset

To delete a user preset,

1. Press **Utility > Imaging Preset Manager**. Select the user preset you want to delete.
2. Press **Delete**. The Delete Preset pop-up menu appears.
3. Confirm that you want to delete this user preset and press **Yes**.

Share My Presets between Versana Premier/Versana Premier LotusSystems

You can share the My Presets you have created between Versana Premier/Versana Premier Lotus systems by exporting/importing the preset(s) you want to share.

To move a user preset from one ultrasound system to another (same software level), first export the user preset(s) you wish to share.

Exporting User Presets

To export a user preset (or presets),

1. Enter **Utility > Imaging Preset Manager**.
2. Insert the media.
3. Press **Export** (on the bottom).
4. An Export Presets pop-up menu appears that indicates:
 - a. Destination location (USB Flash Drive/Hard Disk Drive location).

- b. Preset directory where the preset should be saved (Preset Export).
- c. Available presets on the scanner.

Select the name for the Preset Directory from the Preset Directory pull-down.

5. Select the **User Defined Presets** under **Available presets** on scanner and press **Export**.
6. Upon a successful export, an informational message will pop-up indicating that "1 preset successfully exported." Press **OK**. Then press Exit to close the Export Presets pop-up menu.
7. Press **F3** to eject the media. Take the media to another system and follow the Importing User Presets instructions below.

Importing User Presets

To import a User Preset,

1. Activate the **Imaging Preset Manager** from the **Utility**.
2. Insert the media (Flash Drive, USB Hard Drive).
3. Press **Import**. The Import Presets pop-up appears and displays the Source Directory and Available Imaging Presets.
4. Select the **User Defined Presets** under **Available Imaging Presets** and press **Import**.
If these presets are already on this Versana Premier/Versana Premier Lotus, you will be asked whether you want to:
 - Overwrite this preset (Yes, Yes to All, No or No to All).
 - Rename this preset (Type the new name and press Rename).
 - Cancel.
5. Upon a successful Import, an informational message will pop-up indicating that "1 preset successfully imported." Press **OK**. Then press **Exit** to close the **Import Presets** pop-up menu.
6. Press **F3** to eject the media.

Storing Images and Cineloops

Images and cineloops that are stored during a current examination are displayed as thumbnails on the clipboard.

When an image is stored, all the additional information that is displayed is saved with it (i.e. probe and application selected, image setting, annotations or measurements).

See Dataflow in Connectivity for detailed settings related storing image/Cine.

Image archive is set by the dataflow selected (See Dataflow in Connectivity for more information.)

To print/store an image, **P1** is most commonly used for the primary destination and internal hard drive.

Storing an Image

To store an image:

After The Exam Is Over

1. While scanning, press **Freeze**.
2. Scroll through the cine Loop and select the desired image.
3. Press the appropriate Print key.

The selected image is stored (per your preset instructions) and a thumbnail is displayed on the clipboard.

NOTE

The ultrasound system numbers the images which are saved in the Local Archive (Instance Number). However, Instance Numbers may change or get duplicated when adding/deleting images to the exam. So for identification, the recommendation is to use Content Date/Content Time on the DICOM server instead of Instance Number.

Storing a CINE Loop

A CINE loop is a sequence of images recorded over a certain time frame. The stored CINE Loops are displayed chronologically on the clipboard.

CINE loops can be stored at any time during scanning. You can choose to preview the cine loop before storage and save the cine loop directly, as described below.

The system can be configured to perform either

- Prospective clip: The system begins storing CINE from when you press the Print button, based on the Time Span setting.
- Retrospective clip: The system stores CINE predetermined time before you press the Print button, based on the Time Span setting.

NOTE

The ultrasound system numbers the images which are saved in the Local Archive (Instance Number). But Instance Number may change or get duplicated when adding/deleting images to the exam. So for identification, the recommendation is to use Content Date/Content Time on the DICOM server instead of Instance Number.

Previewing and Storing a CINE Loop

1. While scanning, press **Freeze**.
2. Move the Trackball to activate Cine.
3. Use the trackball or **Frame by Frame** to scroll through the acquisition and find the sequence of interest.
4. Press **Start Frame** or **End Frame** to set the corresponding cineloop boundary to the current frame as necessary.
Rotate **Start Frame** and **End Frame** to trim or expand the cineloop boundaries.
5. Press **Run/Stop** to run the cineloop and then press the print key to store the cineloop.
Cine loops stored on the clipboard are indicated with a movie strip icon.
6. Press **Run/Stop** again to stop the cine loop.
7. Press **Freeze** to return to live scanning.

Depending on whether the system has been configured to enable or disable "Preview Loop before store", the following procedures enable the CINE loop to be stored directly.

Storing a CINE Loop Without Preview

If “Preview Loop before store” is disabled:

1. While scanning, press the appropriate print key.
2. The last valid CINE loop is stored in the archive and a movie clip thumbnail is displayed on the clipboard.
3. Scanning resumes immediately.

Storing a CINE Loop With Preview

If “Preview Loop before store” is enabled:

1. While scanning, press the appropriate print key.
 2. The last valid cine loop is previewed.
 3. Adjust the cine loop, as necessary.
 4. Press the appropriate print key.
- The movie clip thumbnail is displayed on the clipboard.

Data Transfer

The user can select and access the Exam Transfer services from the Exam Data Transfer screen.

- Import
- Export
- Worklist
- Q/R (Query/Retrieve)
- Quick Save

NOTE

Ensure that all patients are exported or backed up BEFORE deleting them.



WARNING

Do not connect or disconnect the probes during data transfer.



WARNING

To avoid data loss, please do not unplug USB during data transfer.

Export/Import

To move exams from one ultrasound system to another system or to back up/retrieve exam information, exam information must be imported/exported.

NOTE

Both database information and images are exported. No data is deleted from the local archive when exporting data.

NOTE

Importing patient records may take more than ten (10) minutes. Please allow sufficient time to import patients.

NOTE

The media used **MUST** be verified **BEFORE** performing Export/Import. Verify the media once each session. If problems are encountered, eject the media and then re-insert the media; then try the Export/Import again.

NOTE

If exporting a previously backed-up exam, the message “Can't Find Source file” displays. The image data has already been removed from the hard disk drive with EZBackup.

NOTE

It is **STRONGLY** recommended to verify files on Eject when using Export.

NOTE

When dealing with a very large dataset (~500GB), the DICOMDIR file may be corrupted by the system during the DICOM Export operation. This happens in extremely rare situations, only with very large datasets, and the patient/exam dataset is not damaged or corrupted; the user may not see the patient/exam listing. In this case, regenerate the DICOMDIR using 3rd party tools. Please contact GE HealthCare service for any assistance.

Exporting Data

To export an exam(s) to a compatible Ultrasound system:

1. Format and label the removable media. Answer Yes/OK to the messages.

NOTE

The system formats the unformatted CD-R/DVD-R automatically when you select Export on the data transfer screen.

2. Press **Patient** and select **Data Transfer**.
3. The Data Transfer screen is displayed. Select Export.
4. “Local Archive-Int.HD” displays on the Transfer From pull-down menu and the patient list included in the Local Archive displays.
5. Select the destination at the Transfer To pull-down menu.
6. Select the patient(s) to export by using the Transfer From search field (the upper field).

Windows commands can be used to select more than one patient.

To select a consecutive list of patients, click the cursor on the first name, move the cursor to the last name, then press and hold down the Shift+right Set key to select all the names.

To select a non-consecutive list of patients, click the cursor at the first name, move the cursor to the next name, then press and hold down the Ctrl+right Set key, move the cursor to the next name, then press and hold down the Ctrl+right Set key again, etc.

Search for patients via the Search key and string.

Or, use Select All to select all patient.

NOTE

Use best judgment when moving patients' images. If there are many images or loops, then only move a few patients at a time.

7. Press **Transfer**. The progress bar displays during the transfer.
8. Press F3 to eject the media. Specify to finalize the media.

NOTE

A dedicated viewer is needed to display exported DICOM or Raw DICOM images on a PC.

Importing Data

To import an exam(s) to another ultrasound system:

1. On the ultrasound system to import to, insert the media.
2. Press **Patient** and select **Data Transfer** (for Traditional Workflow). Press **Patient** and select **Archive** (for Simplified Workflow).
3. The Data Transfer screen displays. Press **Import**.
4. Select the media from the **From** pull-down menu.
5. The **From** search field shows the patients available for import from the removable media just loaded onto the system.
6. Select the patient(s) or the exam(s) from the list to be imported.
7. Press **Transfer**. The progress bar displays during the transfer.
8. Wait for the patient information to be copied to the system. Informational messages appear while the import is taking place.
9. Press **Eject** to eject the media.

NOTE

Use Import to restore EZBacked up and images.

NOTE

Data can be retrieved from the media to the local drive, played back, or exam information can be processed on the system as Raw Data.

NOTE

L6-12-RS Thyroid/Breast Color raw data is not compatible on High and Low configuration. L6-12-RS Thyroid/Breast Color raw data acquired from High configuration cannot be used on Low configuration (will have artifact).

Query/Retrieve (Search and retrieve the data from DICOM device)

NOTE

For Query/Retrieve to find a patient, the patient MUST have a patient ID.

NOTE

Before you retrieve data from the Worklist server, make sure that default IP address is input in the Default Gateway field in **Utility > Connectivity > TCP/IP**.

Query

1. Press **Patient** and select **Data Transfer**. The **Data Transfer** screen displays. (for Traditional Workflow)
Press **Patient** and select **Archive** (for Simplified Workflow).
2. Select Q/R. The patient/exam list in the Local Archive displays in the To section.

NOTE

Only "Local Archive - Int.HD" is enabled for Transfer To.

3. Select the Query/Retrieve server from the From pull-down menu.

NOTE

The server is configured in the Utility screen. Multiple servers are able to be configured.

4. Press **Query** in the transfer **From** section. The Query is performed.
5. The server's patient list displays.

NOTE

Press Query again to refresh the list.

Retrieve

1. Select the patient(s) or the exam(s) to be retrieved from the patient list.
2. Select **Transfer**. Retrieve the data from the Query/Retrieve Server. The progress bar displays during the transfer.

Worklist (Search and Retrieve the Patient/Exam Information)

NOTE

Before retrieving data from the Worklist server, ensure the default IP address is input in the Default Gateway field in **Utility > Connection Manager > Network > Wireless**.

NOTE

Select the patient prior to sending images to a PACS device.

NOTE

User can configure the Worklist to the User Defined Key through **Utility > System > User Configurable Key > User Defined Key** (Traditional workflow only).

1. Press **Patient** and select **Data Transfer**. The Data Transfer screen displays.
2. Select **Worklist**. The patient/exam list in the Local Archive displays.

NOTE

Only “Local Archive - Int.HD” is enabled for Transfer To.

3. The previous Worklist used is displayed on the monitor display. Press **Query** to refresh the list.

NOTE

The worklist server is configured on the **Utility** screen. Multiple servers can be configured.

The auto-refresh worklist can be enabled/disabled on the **Utility** screen. The system automatically refreshes the list when the exam data transfer accesses the Worklist server or changes the Worklist server.

4. Select the patient(s) or the exam(s) from the list.

NOTE

User can also double click to enter a selected patient or exam.

5. Press **Transfer**. The progress bar displays during the transfer.

USB Quick Save


Quick Save easily sends Images to a USB flash drive storage.

The Images are stored either in JPG or WMV format.

USB Quick Save Setup

In **Utility** > **Connectivity** > **Service**, select **USB Quick Save**.


Figure 3-6 USB Quick Save



After The Exam Is Over

The USB Quick Save service can also be assigned to the **Print** keys via the **Utility > Connectivity > Button** preset menu.

Figure 3-7 Assigning USB Quick Save to Print Keys



BluetoothSave

BluetoothSave can easily send images to a Bluetooth storage device (Android devices only). User can press the configured hotkey to send current image or cineloop to folders in the storage device.

Bluetooth Device Configuration

To add a new Bluetooth device,

1. Connect bluetooth adapter to the system and reboot the system.
2. Press **Add New Device** in **Utility > Connectivity > Bluetooth**.
3. Follow the steps in the dialogue to complete configuration.
4. Type the device name in the **Name** field. Android device name will displayed as default name.

Table 3-1 Device

Preset Parameter	Description
Add New Device/Remove	Press Add new device to add a new device; press Remove to delete a device.
Properties: Name	Type the name of the device.

Preset Parameter	Description
Properties: MAC Address	Unique network card address. NOTE Only available for MyComputer.
Properties: Date Added	Display the new device added date.
Default Bluetooth Device	Select default Bluetooth device from the list.

System hotkey configuration

1. Assign **BluetoothSave** as **Numeric Hot Key** or **Alphabet Hot Key** though **Utility > System > User Configurable Key**.
2. Return to scan mode and press **Freeze** key, then press **BluetoothSave** hotkey.
3. Press **OK** to complete.

Send To (Send the image to the DICOM Device)

Send To sends the selected exam for a patient interactively to a destination DICOM device configured on the system. An exam in this case includes its images and any corresponding Structured Report.


1. Search and select the patient and press **Review**. The Exam View screen displays.
2. Select the exam which has the images and press **Send To** icon.

NOTE

You can only select "Local Archive - Int.HD" for Workflow.

3. The "Send To" Dialogue box displays.

Figure 3-8 Send To Dialogue



Select the destination device and press **OK**.

NOTE

The destination device is configured in the Utility screen. Multiple devices are able to be configured.

4. The successful/unsuccessful message is displayed at the top of the screen.

Selectable Send To

Selectable Send To sends the selected images for an exam and to interactively Send-To the exam from Local Archive or DICOM Read to a destination DICOM device configured on the system.

1. Select a patient or an exam in the **Archive** menu.
2. Go to **Active Images** (for Traditional Workflow).
3. Select an image (images) and press **Send To**. Send To Dialogue displays.

NOTE

If the image is not selected, warning dialog is displayed and no image is sent.

4. Select a destination from pull-down menu and press **OK**.

"**Selected image(s) is (are) send to**" message displays on the status bar.

Using the DICOM Spooler

To monitor/control DICOM jobs, press F4. You can view, resend, redirect, and delete images from the DICOM spooler by selecting a job, then specifying the action to be performed on this job.

NOTE

If you find a failed job(s) in the Spooler, please remove the failed job(s) from the Spooler.

Table 3-2 Spooler status description

Status	Description
Hold	Waiting for user activity. Select Resend or Send-To to complete the job.
Pending	Waiting for the previous job(s) to finish (a previous job may be Active or Pending). No user interaction required.
Append	Not completed. Example 1: Direct Storing job. Waiting for more images or the end of the exam (by selecting New Patient or End Current Patient). Example 2: Print job with 3x3 images has only 8 images. Waiting for one more image or the end of the exam (by selecting New Patient or End Current Patient).
Active	Signifies network activity (or connection attempt).
Success	Sent successfully.
Failed	Unsuccessful job attempt. Job stays in spooler. Select Retry or Delete to complete the job.
Done	Finished successfully.

Image/Data Management

Reviewing Patient Images

Retrieving and editing archived information contains below section:

- Searching for a patient
- Reviewing a patient exam
- Reviewing an image
- Deleting a patient, exam, or image

Clipboard

The clipboard displays thumbnail images of the acquired data for the current exam. Images from other exams are not displayed on the current patient's clipboard.

All of the images can be viewed in the Active Images screen, available from the display or from the Archive menu.

Saving the image /cine to the Clipboard

The active image/cine is stored and placed on the clipboard when you press the Store key. The clipboard contains preview images with enough resolution to clearly indicate the contents of the image. CINE Loops are indicated by a movie clip icon.

Previewing Clipboard Images

1. Select the **Cursor** key to obtain a cursor arrow.
2. Move the **Trackball** to position the pointer over the clipboard image you want to recall.
3. An enlarged preview of the image is displayed on the left-hand side of the monitor.

Recalling Images from the Clipboard

To recall images from the clipboard,

1. Select the **Cursor** key to obtain a cursor arrow.
2. Move the **Trackball** to position the pointer over the clipboard image you want to recall.
3. Press **Set** to recall the image.

To delete an image from the clipboard

1. Select the **Cursor** key to obtain a cursor arrow.
2. Place the cursor on the clipboard image you want to delete, then press **Set** to select the image.
3. Place the cursor on the **Delete** icon and press **Set**.
A warning message is displayed asking the user to confirm the action to perform.
4. Select **OK**.

SaveAs

Images and cineloops can be saved to a removable media storage to View on a Windows PC in the following standard formats:

- Still images: JPEG, DICOM and RawDICOM (Raw data + DICOM)
- Cineloops: AVI, DICOM and RawDICOM (Raw data + DICOM)

Using Save As

To save images to media:

1. Insert the media into the drive or connect the USB drive to the system.

NOTE

If the media is not formatted, it will be formatted when you select **Save As**.

2. On the scan screen, press the **Cursor** key. The arrow cursor displays.
3. Place the cursor on the image or CINE Loop in the clipboard to be saved and press **Set**. The image displays on the screen.
4. Select **SaveAs** in the lower of the screen. The Save As menu appears.

NOTE

If the image is saved as an .avi file, run the CINE Loop before you select SaveAs.

NOTE

A 2D cine loop image cannot be saved as a .jpeg file.

Figure 3-9 Save As Menu

5. Select the media from the **Save in archive** pull-down menu.
6. Folder name: Create the folder for the saved file.
 - Default is blank (The folder is not created)

NOTE

The folder name can not be edited when the folder is opened.

7. File Name: The name of the file is automatically completed, but a desired file name can be typed as well.

NOTE

DO NOT use the following special characters when saving images: !, @, #, \$, %, ^, &, *, (,), |, :, ;, <, >, ?, /, ~, [,], {, }, and Yen sign.

8. Store: Select Image only or Secondary capture.
 - Image only: Saves only the ultrasound image area
 - Secondary capture: Saves the ultrasound image area, title bar, and scan information area. Not available for DICOM or Raw DICOM images.

NOTE

If “WMV” is selected for the Save as type, Secondary capture is disabled.

9. Compression: Specify Compression.

- None
- Rle
- Jpeg

NOTE

If “WMV” is selected for the Save as type, Compression is disabled.

10. Quality: Specify image quality (between 10-100). A high quality setting gives a lower compression.

NOTE

If “WMV” is selected for the Save as type, Quality is disabled.

11. Save as type: Select one of the following.

- Raw DICOM: saves the still image or CINE Loop in both GE raw format and DICOM format.
- DICOM: saves the still image or CINE Loop in pure DICOM format.

NOTE

DICOM Images shall be saved as Jpeg image. DICOM loops shall be saved as AVI or WMV format.

- AVI: Saves the CINE Loop in AVI format.

NOTE

Store “Image Only” is available if you select AVI for Type.

- JPEG: Saves the still image in jpeg format.
- WMV: Saves the CINE Loop in wmv format.

NOTE

WMV type is only available with CINE loop image.

NOTE

The Save button is disabled when you select “AllFiles.” Select each Save as type when saving the data.

If you want to see all data saved onto the SSD, select “AllFiles(.*)”. All the data names display in the window.

12. For images transferring to USB, press **Save**.

The images are saved directly to the USB drives storage whenever you press Save.

If saving to a CD/DVD, select “For Transfer to CD/DVD,” the image is saved to the local drive buffer.

- If there is not enough free space in the destination to save all selected images, a warning dialog will appear.
- If the same file name exists in the destination, the warning dialog displays.

OK: Overwrite file and continue to save selected images.

Cancel: Cancel.



WARNING

Buffer images for transfer to CD/DVD. Images are not cleared out of the image buffer after you start a new exam or new patient. It is the responsibility of the user to clear the image buffer before storing new patient data so that the CD/DVD does not contain images from multiple patients.

13. Repeat steps 4-12 for each additional image to be stored to the CD/DVD.
14. After adding all of the images/loops to save and ready to write to the media, transfer all the images at the same time. Press **Save As > Transfer To CD/DVD**.

A progress bar displays, notifying that the “Media transfer is in progress.”

If the total transfer size is bigger than CD/DVD free space size, then only the files that can be copied to CD/DVD are transferred. After the copy is finished, a warning dialog displays showing the total required file size and transferred file size. Press **OK**, insert a new CD/DVD, and press **Transfer to CD/DVD** again.

15. To not save the image(s) to CD/DVD, select “**Delete Files for Transfer**” All images are deleted.
16. Press **Eject** to eject the media. Select CD/DVD Recordable or USB drive. Select **Yes and Verify files** for CD/DVD. This compares the expected number of files with the actual number of files on the media. The files are also checked to ensure that they are readable.

NOTE

The Report Save As feature is somewhat different. As soon a report is selected to be saved, the report is saved.

NOTE

If a 3D image is saved as an WMV file, an annotation text “COMP” appears at the top of the saved image which represents the compressed image.

NOTE

Time line image can be saved as multi frames image with SaveAs.

Table 3-3 Save As Formats

	.wmv format	MPEG Vue in Data Transfer
B, B+CF	Multi frames	Multi frames
B mode + Doppler mode	Multi frames	Single frame
B+M	Multi frames	Single frame
3D	N/A	Single frame

After The Exam Is Over

Verify the saved image works correctly on the Windows PC. If the image does not work, please save it again on the ultrasound system.

'SaveAs' Images - Traditional Workflow only

Multiple images can be selected to save at one time by selecting **'SaveAs' Images (Send To for Simplified workflow)** in the **Active Images (Image Review for Simplified workflow)** screen.

Features are almost the same as the **Save As** feature. See **SaveAs** on page 172 for more information.

NOTE

It is suggested to save images page by page. Saving several images or raw data at once can take a longer time.

NOTE


If the image has a filmstrip icon, this indicates a CINE Loop, which can be saved as a .wmv or .mp4 file; single images can be saved as a jpeg file.

NOTE

'SaveAs' Images function doesn't support images which are query/retrieved.

1. In the Active Images screen, place the cursor on the image or CINE Loop to be saved and press **Set**. Multiple images with multiple pages can be selected.
2. Press **'SaveAs' Images (Send To for Simplified workflow)** on the monitor display or the Touch Panel. The Save As (Send To for Simplified workflow) menu appears.

Figure 3-10 Save As menu



3. Ensure that **Jpeg & WMV** or **Jpeg & MP4 (JPG / AVI, JPG / WMV, DICOM or JPG / MP4 for Simplified workflow)** is selected, then press **Save (Send for Simplified workflow)**.

NOTE

If you are saving to USB, images are transferred as soon as you press Save; if you are saving via Transfer to CD/DVD, you need to save images to the hard drive, then images are transferred when you select Transfer to CD/DVD. See below.

Storing Images with More Resolution

To store images with more resolution than is available with the JPEG selection, select Save As and select AVI as the Save As Type. Single images can be saved as .avi files.

Table 3-4 Store Options

Image Type	Store as Image Only	Store as Secondary Capture
CINE Loop	Provides a loop of just the image (no title bar and scan information).	Provides a single image of the video area. DO NOT USE, UNCLEAR WHICH IMAGE FROM THE LOOP IS SAVED.
Still Image	Provides a single image (no title bar and scan information).	Provides a single image of the video area.

Data Backup

The backup and restore procedures described in this section are divided into two parts. The first part describes procedures to backup and restore patient data. The second part describes procedures to backup and restore system and user-defined configurations.

The Backup/Restore function enables the user to:

- Copy/Restore the patient archive.
- Copy/Restore the system configuration. The Copy/Restore system configuration feature enables the user to configure several units with identical presets, providing that the units have the same software version.

Depending on the system, you can use either a CD-R, DVD-R, USB Flash Drive, or USB Hard Disk for system backup/restore. For the sake of simplicity, we have used the CD-R in the following examples.

NOTE

The system ONLY supports CD-R / DVD-R and DOES NOT support CD-RW / DVD+R.

**WARNING**

GE is not responsible for lost data if the suggested backup procedures are not followed and will not aid in the recovery of lost data.

**WARNING**

The Versana Premier/Versana Premier Lotus is not intended to be used as a storage device; backup of the Patient and Image Database is your institution's responsibility. GE is NOT responsible for any lost patient information or for lost images.

**WARNING**

The system crash can cause the internal SSD corruption. The SSD is not considered a permanent storage device. Backup data on a regular basis.

**CAUTION**

To minimize accidental loss of data, perform EZBackup and Backup on a regular basis.

1. First, perform EZBackup to save the images.
2. Next, perform Backup at Utility -> System -> Backup/Restore. Enable the following checkboxes under Backup:

- Patient Archive
- User defined configuration
- Service

**CAUTION**

Archived data is managed at the individual sites. Performing data backup (to any device) is recommended.

**CAUTION**

Make sure to verify the media after writing of data, such as EZBackup, SaveAs or Export. Verifying media requires additional time, which varies depending on the amount of data backed up or exported.

**CAUTION**

Before deleting a patient or image from the patient screen, make sure you have saved the data by EZBackup/Backup or Export and verify that the media transfer of data was successful.

EZBackup and EZMove

EZBackup or EZMove allows you to manage hard disk space (move images off the hard drive) while maintaining the patient database on the scanner, as well as to back up the patient database and images.

**HINTS****PLEASE READ THIS**

Ensure that you have established a data management protocol for your office/institution. You **MUST** manage the backup media by keeping a log and by creating a media filing system.

For example, if you need to back up 500 MB/day, or 2.5 GB/week, then you need to back up 5 CDs/week, or ~250CDs/year.

Generally speaking, you should back up the system when you have 10 GB of images to back up.

You should assign the person who is in charge of performing the backups. Backups will vary by the volume of your work. You need to track how long it takes your office/institution to get to 10 GB, and set the back-up parameters accordingly.

Your office/institution needs to determine your backup strategy, for instance, backup weekly and move monthly. It should be an easy strategy to perform and to remember. And follow this same strategy/schedule consistently.

It's also useful to keep your more recent information on the hard drive since it's easier to recall that way.

**CAUTION**

You can still do a backup/move daily; but **ALWAYS** do a patient archive backup after each move.

**CAUTION**

Only cancel the backup/move in case of an emergency. The system completes backing up the current media and then cancels the operation.

**CAUTION**

When EZBackup requires more than one disk (CD-R or DVD-R) for backup, a message appears when the first disk is full. If you select "Cancel" to stop the backup procedure and later try EZBackup again, all the data may not be backed up.

Select "Full Backup" on the first EZBackup wizard screen if the last time you were performing EZBackup you selected "Cancel".

**CAUTION**

If you use EZBackup or EZMove as a "true" patient archive, you must maintain a separate backup of the patient database (Patient Archive and Report Archive). If for any reason the Local Archive - Int HD gets corrupted or the base system software has to be reloaded, then the patient archive is the **ONLY** way to rebuild the EZBackup and EZMove patient archive.

**CAUTION**

DO NOT turn off the power while EZBackup is running. The data may be lost. It may take several hours for EZBackup to finish, depending on the amount of data being backed up.

The following may give the impression of a lockup, but EZBackup is continuing in the background.

- The progress bar does not move.
- The screen may become white.
- The hourglass icon keep turning.

**CAUTION**

NEVER restore the patient archive from media made previous to the last move.

NOTE

EZBackup/EZMove saves data as RAW data. If you import data to the system, you can modify the image data.

NOTE

To display exported Raw DICOM images on a PC, you need the dedicated viewer.

NOTE

When backing up or moving reports, do not use EZBackup or EZMove, EZBackup or EZMove cannot back up or move reports.

NOTE

"Archived" information is saved to each exam during EZBackup. When you perform EZBackup, the system backs up the exams except for the archived exam.

NOTE

EZBackup/EZMove cannot span a single image across two (2) or more media. Therefore, if EZBackup/EZMove encounters an image that is greater than the capacity of the media, it skips the oversized image.

NOTE

EZBackup/EZMove does not store images to media in sequential order. Instead it maximizes the most amount of images per media.

NOTE

If the system locks up during the media auto format process, shutdown the system by holding down the power button and reboot the system. After the system is up, replace the media to a new one and execute EZBackup or EZMove again. To avoid trouble such as data loss, do not reuse the failed media for any other function.

NOTE

If you try exporting a previously backed-up exam, the message "Can't Find Source file" displays. The image data had already been removed from the hard disk drive with EZBackup/EZMove.

Basically, when you perform the EZBackup or EZMove procedure, you insert the media (or connect USB HDD if applicable), the system backs up/moves the images, and creates a reference between the patient database and the media's volume.

EZBackup/EZMove can take up to 20 minutes (or longer, depending on the size of the backup). Make sure to schedule this at the same time daily, when no patients are scheduled.

1. Prepare unformatted media or the USB HDD before starting EZBackup/EZMove.
2. Specify the EZBackup/EZMove setup on the **Utility > System > Backup/Restore** page.
3. To start the EZBackup/EZMove procedure, go to the **Patient** menu and select EZBackup/EZMove. The EZBackup/EZMove Wizard starts.

NOTE

If you use the USB HDD, some wizards and the pop-up messages DO NOT appear.

4. Verify the information on the first page of the EZBackup/EZMove Wizard, then press Next.

Full backup options display on the first page of the EZBackup wizard. If you want to backup all of the exams in the range (even if the exam was previously backed up, check this option). If you uncheck this option, the system only backs up exams which have not yet been backed up.

EZBackup does not back up the exams which were previously backed up once by EZBackup or Export.


NOTE

You can set the range in **Utility > System > Backup/Restore > Move files older than in days for EZMove**.

NOTE

If you update an exam which is already backed up, the exam is also backed up.

Figure 3-11 EZBackup/EZMove Wizard 1



5. Verify the information on the EZBackup/EZMove Wizard, Page 2. The backup may span multiple media. This page tells you how many media you need to do this backup. After

After The Exam Is Over

you have gathered the media (allow for one extra media, just in case), you are ready to begin the backup. Press Next.

Free Space/Total Size: tells you the size of the data you have selected to store/and the total size of the USB Hard Drive storage media. If the storage capacity of the USB HD is insufficient, you will see the message, "Selected Location does not have enough free space."


NOTE

The calculation for the number of backup CD is only an estimate. Allow for one additional CD when performing an EZBackup/EZMove.

NOTE


This message appears if you press Next without inserting the backup media: "Please insert a blank media...". Insert the media and continue.

Figure 3-12 EZBackup/EZMove Wizard 2



6. A pop-up message appears that provides you with the media label. Label the media, then insert the media. Press OK.

Figure 3-13 Insert Media Message




- a. Ensure that you label the media with not only the volume name indicated on the Insert Media Message, but with the name of the ultrasound system where this backup/move procedure was done.
- b. Update the EZBackup/EZMove log with this information the volume information and the location of the media.
- c. After the backup/move has been completed, file the media.

Table 3-5 Typical EZBackup/EZMove Log

Date	Scanner ID Name	Backup Images Y/N	Older than _ Days	Move Images Y/N	Media Label (and Scanner ID)

7. The system will check whether the media is inside the drive.
 - If the media is in the drive, the window automatically disappears when the checking is completed. The EZBackup progress then begins.
 - If the media is not in the drive, a pop-up message appears notifying that there is no media in the drive. Press **OK**.


Figure 3-14 Insert Disc Error



A pop-up message appears that provides you with the media label. Label the media, then insert the media. Press **OK**.

8. The status menu appears. When the backup/move has been completed, press **Next**.

Figure 3-15 EZBackup Wizard 3



NOTE

When/if you need to insert the next media, a message appears providing you with the media label. Label the media, then insert the next media and press **OK**.

9. When the backup is complete, the completed wizard page appears. Press **Finish**.
10. Do a patient archive after each EZBackup/EZMove.

We recommend attaching the patient list to the EZBackup/EZMove media.

NOTE

Use Import to restore EZBackup images.

To Review EZBackup/EZMove and Export Images

You can review backed up media via the Patient Menu, Import, and the DICOM Read dataflow.

If you review EZMoved image,

1. Select the patient on the Patient Menu (on the same system where the EZMove was performed).
2. Insert the media volume indicated on the Patient Menu
3. View the exam from the media.

NOTE

You may need to insert a media volume prior to or after the recommended media.

NOTE

If the patient is split over multiple media, images on the previous or next media are displayed as triangles.

NOTE

To view the whole patient on the system, use Import, from as many media as you have for that patient. However, take care not to import studies over existing studies; duplicate or missing images may result. Delete the existing exam first.

Backup and Restore

To minimize accidental loss of data, perform backup of the patient archives stored on the local hard drive DAILY as described in this section. Use a formatted Backup/Restore disk to back up patient archives from the hard drive, using the backup procedure described in this section. Data from the Backup/Restore disk may be restored to the local hard drive using the restore procedure.

NOTE

To perform backup and restore procedures, you must login with administrator privileges.

Backup procedure

Back up patient data AFTER you've archived (via EZBackup/EZMove) images so that the pointers to the patient's images reflect that the images have been moved to removable media and are no longer on the hard drive.

1. Insert a media into the drive or USB device into a USB port.
2. On the **Archive** screen under dataflow, select **Local Archive - Int. HD**.
3. Go to **Utility > System > Backup/Restore**.
4. In the Backup list,
 - Select **Patient Archive** to backup the patient records.
 - Select **User Defined Configuration** to copy system settings and user presets.
 - Select **Service** to backup service presets.

NOTE

The detailed section of this menu decouples the user defined configuration above. This allows you to selectively restore what you want to restore across multiple machines.

5. Specify where to save data in the media field.
6. Select **Backup**.
The system performs the backup. As it proceeds, status information is displayed on the Backup/Restore screen.
7. At the end of the process, the Backup completed message is displayed on the monitor.
Press **Eject (F3)** for eject media/disconnect USB device.
8. Make sure to physically label the media. An identification of the system should also be noted on the media and a backup log should be kept.

File the media in a safe place.

Restore procedure



CAUTION

The restore procedure overwrites the existing database on the local hard drive. Make sure to insert the correct media.

You cannot restore the data between systems with different software versions.



CAUTION

To avoid the risk of overwriting the local patient and report archives, DO NOT check Patient Archive when restoring user-defined configurations.

1. Go to **Utility > System > Backup/Restore**.
2. In the Restore list,
 - Select **Patient Archive** to backup the patient records.
 - Select **User Defined Configuration** to copy system settings and user presets.
 - Select **Service** to backup service presets.
3. In the Media field, select the appropriate Source device.
4. Select **Restore**.

The system performs the restore. As it proceeds, status information is displayed on the Backup/Restore screen.
5. System restarts automatically when restore is done.

Backup and restore strategy: user-defined configurations

In addition to generating a safety copy, the backup/restore function of the user-defined configuration (presets) can be used to configure several Versana Premier/Versana Premier Lotus systems with identical presets (preset synchronization).

Preset synchronization

The procedure for preset synchronization of several scanners:

1. Make a backup of the user-defined configurations on a removable media from a fully configured Versana Premier/Versana Premier Lotus system.
2. Restore user-defined configurations from the removable media to another Versana Premier/Versana Premier Lotus system (you can restore all the user-defined presets or select specific presets to restore via Detailed Restore).

Configuring Connectivity

You use Connectivity functionality to set up the connection and communication protocols for the ultrasound system. This section gives an overview of each of the Connectivity functions.

Connectivity Functions

To set up your institution's connectivity, you must login with administrator privileges.

1. **TCP/IP:** allows you to configure the Internet Protocol including IP settings and wireless configuration.
2. **Device:** allows you to set up devices.
3. **Service:** allows you to configure a service (for example, DICOM services such as printers, worklist, and other services such as standard print) from the list of supported services. This means that the user can configure a device with the DICOM service(s) that particular device supports.
4. **Dataflow:** allows you to adjust the settings of the selected dataflow and associated services. Selecting a dataflow customizes the ultrasound system to work according to the services associated with the selected dataflow.
5. **Button:** allows you to assign a pre-configured output service (or a set of output services) to the Print keys on the control panel.
6. **Removable Media:** enables formatting (DICOM, database, or blank formatting) and DICOM verification of removable media.
7. **Miscellaneous:** allows you to set up the patient/exam menu options, print and store options, patient/exam Message options, Other ID Options and columns in the examination listing.
8. **Bluetooth:** allows you to set up Bluetooth device connection.

Configure these screens from left to right, starting with the **TCP/IP** tab first.

NOTE

The ultrasound system is pre-configured for many services, with default settings selected. You can change these services and settings as needed.



CAUTION

You must restart the ultrasound system after making any changes to connectivity settings in the Utility menus. This includes any changes on the **TCP/IP** or **Dataflow** setup screens.


Imaging Insights

Imaging Insights automatically collects DICOM data from GE HealthCare and other vendors' ultrasound equipment or equipment fleet and displays system utilization and operator usage insights in a plotted dashboard. Operational insights include exam volumes, first and last exam time, probes utilization, exam type; operator usage data includes length of exam, scan mode, probes, and exam type. It helps optimize system and probe fleet investment plans, identify staff assignment and training needs, and monitor variability of staff usage patterns.

After The Exam Is Over


1. Press **Utility > Connectivity > TCP/IP > Imaging Insights** to enter **Imaging Insights** configuration window.

Figure 3-16 Imaging Insights button



2. Configure the server information in the pop-up window and then check **Send data to Imaging Insights** check box.
3. Click **Save Settings** to save the configuration and send data to **Imaging Insights** server.
 - a. If the network connection to the server passes, data will be sent to server successfully.
 - b. If the network connection to the server fails, a question will pop out for user.

Figure 3-17 Imaging Insights question



Typical failed causes:

- Network cable not connected.
- Configuration error(s).

NOTE

If the user is doing the exam with patient information created, the data will be sent to the server immediately when the user clicks **End Exam**. If the user is scanning without creating patient information, the data will be sent to the server every two hours before system shutdown.

NOTE

User can press **View Statistics** button to check whether the data is sent successfully or not.

Anti-Virus Software Note

Anti-virus software IS NOT present on the Versana Premier/Versana Premier Lotus system. Since the ultrasound is already protected against attack by the measures listed below, no Anti-virus software is deemed necessary.

- Only communication ports required for system operation are enabled.
- Only operating system services required by system application software are enabled.
- Software programs CANNOT be loaded onto the ultrasound system (e.g., e-mail, web browser, etc.).
- An auto-executable file CANNOT be run automatically on the ultrasound system.
- The ultrasound system software includes the latest MS Windows security protection.

We have worked diligently to develop a combination of the safety measures noted above and the security standards of Windows 10 Enterprise 2016 LTSC to provide a degree of safety against Viruses, Worms, Trojan Horses, etc., especially for a system used in a professional hospital grade networking environment that also typically features its own sufficient safety measures.

Printing Options

Setting up Digital Peripherals

You set up digital peripherals from the **Utility > System > Peripherals** menu.

NOTE

Printing using a standard printing service overrides the orientation and N-up feature of the printer preferences. Printer preferences are set up in the printer folder (via **Utility > System > Peripherals**. Select **Properties** under **Standard Printer Properties**).

Digital Printer Setup

There are two steps to do when setting up a digital printer: follow the procedure below for each printer, then set up specific properties for each printer if you need.


Follow this procedure for each printer:

1. Select **Utility > Connectivity > Service**. Add the **Standard Print**. Highlight **Standard Print** in the **Service** list. Select the printer from the **Printer** pull-down **Properties** menu. For the UP-D898 printer, select "Portrait" as orientation.
2. Type the printer name in the **Name** field. This name is used on the **Button** screen. After you select the printer from the Printer pull-down Properties menu again, it turns white. Press **Save**.
3. Select **Button**. Select the appropriate print key (**Print1**, **Print2**, **Print3**) from the **Physical Print Buttons** section as shown in the figure below. Select the printer from

After The Exam Is Over

the MyComputer column and press >> to move it to the Printflow View column. Press **Save**.


Figure 3-18 Digital Printer Setup



Sony UP-D898 Instructions

Follow these steps to set up the paper size of Sony UP-D898 printer.

1. Press **Utility > System > Peripherals**. Select the UP-D898 from the pull-down menu under **Standard Printer Properties**. Click **Properties**.
2. Click **Printer > Printing Preferences** on the menu of Properties Window.



3. Select Paper Size. Press **Apply**. Press **OK**.
4. Press **Save**, then **Exit**.

Setting up the Off-line paper printer

You can connect an off-line paper printer via the USB connection.


**CAUTION**

ONLY plug in devices to the USB ports located at the rear of the system WHILE the ultrasound system is NOT powered up. If you plug in a device while the ultrasound system is powered on, your system may become unusable.

**CAUTION**

DO NOT place an off-line paper printer inside the patient environment. This assures compliance to leakage current.

Figure 3-19 Patient Environment

**NOTE**

The printer driver is customized for the ultrasound system at the factory; you do not need to change the settings.

1. Connect the printer to the USB port.
2. Go to **Utility > Connectivity > Service**. Select the printer and press **Add** to add the printer.
3. Select the printer from the Printer pull-down Properties menu.

NOTE

After selecting the printer, the field turns white.

4. Set the following parameters in Properties:
 - Rows=3

After The Exam Is Over

- Columns=2
- Orientation=Portrait
- Right Margin (mm)=10

5. Type the printer name in the Name field.

NOTE

This name is used on the Button screen.

6. Press **Save**, then select the Button tab.
7. Select the **Print** key from the control panel.
8. Select the printer from **MyComputer** column and press >> to move it to the **Printflow View** column.
9. If you want to assign this printer to the Standard Print Button on the Active Image Screen, select this printer at the Active Image Printer section.
10. Press **Save**.

NOTE

If you want settings other than 1 image per sheet, or 2x3 sheet, or to improve the image quality, refer to the manual that came with the printer.

Setting up the Printer to Print Reports

To set up the Off-Line Printer to print reports,

1. Enter **Utility > System > Peripherals**.
2. Select the printer from **Default Printer** pull-down menu.
3. Press **Save**.
4. Press **Print** on the **Report** screen to print the report.

Electronic Documentation (Not applicable for China)

Documentation Distribution

The manuals for user will be supplied via several types: hardcopy, softcopy in media and softcopy inside system software.

- Hardcopy
 - User Manual (translated)
 - Advanced Reference Manual (English and French)
 - Release Notes and Workarounds (translated)
 - Basic Service Manual (English only)
 - AIUM Acoustic Output Booklet (USA only, hard copy)
- Electronic media. You can view user documentation (all languages) on a PC or on the Ultrasound Scanner via the Customer Documentation media, which includes:
 - User Manual (translated)
 - Advanced Reference Manual (English and French)
 - Release Notes and Workarounds (translated)
 - Basic Service Manual (English only)

Using Online Help Via Help

Online Help is available via the **Help** key. After pressing **Help** key, Help screen displays the User Manual in PDF format. The PDF reader allows to browse, search and zoom into the page of User Manual.

Viewing Online Help in a Language Different from the System Language

Go to **Utility > System > General > Manual Language**, select the preferred language to view Online Help.

If the translated Online Help is not available, you will be directed to select a different language.

Translated Online Help files can be updated via the eIFU USB Flash Drive. Please contact your Applications/Field Service Representative to order an updated eIFU Kit.

Go to **Utility > System > General > Manual Language**, select the preferred language to view Online Help.

If the translated Online Help is not available, you will be directed to select a different language.

Translated Online Help files can be updated via the eIFU USB Flash Drive. Please contact your Applications/Field Service Representative to order an updated eIFU Kit.

Browsing through the Bookmarks

Online Help is organized with individual chapters, sections, and pages. Click on the plus (+) sign next to the chapter title to open up that chapter. Click on the section you want to view to open up that section.

Searching for a Topic

To search for a specific topic, click on the Search tab in the tool bar of the Help page. Type in the topic name to search.

Figure 3-20 Search for a Topic



Zoom into the Page

To zoom into a specific page, click on the Zoom icon in the tool bar of the Help page.

Figure 3-21 Zoom into the Page



Using the Index

You can look for topics by using the Index. Press the Index tab, then use the scroll bar to look up a topic.

Exiting Online Help

To exit Online Help, press the 'X' in the upper, right-hand corner of the Online Help window.

Electronic media

Accessing Documentation Via a PC

To view user documentation on a PC,

1. Plug in the media into the media drive.
2. Open the media drive.
3. Double click on gedocumentation.html.
4. Click the item in desired language to view the user documentation.

To close the window, click on the 'X' in the upper, right-hand corner of the browser window.

NOTE

If your PC does not have Adobe Reader, a free download is available on the Adobe web site at <http://www.adobe.com>.

Updating Documentation on the Ultrasound Scanner Via the media

The latest version of the User Manual is available to the user on the USB Flash Drive.

In order to update the system User Manual to the latest version,

1. Power down the ultrasound system and insert the eIFU USB Flash Drive into a rear USB port.

NOTE

Ensure that the system is USB Device Enabled (check setting on System Admin Utility Page).

2. Power on the ultrasound system. Online Help files will be loaded onto the ultrasound system automatically, following several screen prompts:
 - a. Select **Install SW...** on the Start Application screen.
 - b. Press **OK** on the StartLoader screen to continue.
 - c. Press any other key to continue when the first C:\windows\system32\cmd.exe screen displays.
 - d. When a C:\windows\system32\cmd.exe window pops up indicating the installation complete, press any key to restart the system and remove the eIFU USB Flash Drive from the system.

System Care and Maintenance

The user must ensure that safety inspections are performed at least every 12 months according to the requirements of the patient safety standard IEC 60601-1. Refer to the Chapter 10 of Service manual.

Only Service Representatives are allowed to perform the safety inspections mentioned above.

Technical descriptions are available on request.

To ensure that the unit constantly operates at maximum efficiency we recommend that the following procedures be observed as part of the customer's internal routine maintenance program.

Contact the local Service Representative for parts or periodic maintenance inspections.

Expected Service Life Description

The expected service life for the Versana Premier/Versana Premier Lotus system and probes is identified in this table:

Table 3-6 Expected Service Life

Equipment / Accessory	Expected Service Life
Versana Premier/Versana Premier Lotus system	The expected service life for the Versana Premier/Versana Premier Lotus is seven (7) years from the manufacturing date under the provision of regular maintenance by authorized service personnel.
Probes	The expected service life for the probes meets or exceeds five (5) years from the date the probe is placed in service, under the provision that the customer follows the care instructions provided on the Probe Care Card / Accompanying Versana Premier/Versana Premier Lotus Instructions for Use.

Maintenance Schedule

Follow this Maintenance Schedule to maintain optimum system function and patient care:

Table 3-7 Maintenance Schedule

Monthly	Weekly	Daily	After Each Patient
Inspect the following on a monthly basis: - Connectors on cables for any mechanical defects - Entire length of electrical and power cables for cuts or abrasions. - Equipment for loose or missing hardware. - Control panel and keyboard for defects. - Casters for proper movement and locking operation.	Clean and disinfect the following on a weekly basis: - Console - System Cabinet - Removable Trackball/Trackball - Air Filters (weekly, or as needed) - Footswitch - B/W Printer - Monitor - Operator Controls panel	Clean and Disinfect the following areas where Cross Contamination can occur: - Operator Panel - Monitor - Front and Rear Handles	Clean and Disinfect the following after each patient: - Probe - Biopsy Bracket, as applicable. Additionally, Clean and Disinfect any area on the system that has visible contamination from the previous exam NOTE Biopsy Accessories must be cleaned and disinfected or disposed of after each patient. Refer to the Probes Chapter, for instructions on probe cleaning, disinfection and inspection.

**CAUTION**

To avoid electrical shock hazard, do not remove panels or covers from console. Servicing must be performed by GE HealthCare qualified service personnel. Failure to do so could cause serious injury.

If any defects are observed or malfunctions occur, do not operate the equipment but inform a qualified service person. Contact a Service Representative for information.

**WARNING**

Do not perform system care and maintenance in the patient environment.

NOTE

Frequency of the cleaning is dependent on the environment.

Failure to perform required maintenance may result in unnecessary service calls.



CAUTION

When performing cleaning and disinfection procedures, to prevent the risk of system damage, always observe the following precautions:

- Use only cleaning and disinfection materials and solutions as recommended in the procedures described below.
- Never use thinner, benzene, ethanol or methanol alcohol, abrasive cleaners, or other strong solvents, as these may cause damage to the cabinet or monitor. Only use isopropyl alcohol, when instructed to do so.
- Do not spray any liquid directly onto the ultrasound system covers, monitor or keyboard.
- Do not allow any chemical agents to drip or seep into device openings or connections.
- Use only recommended cleaners or disinfectants on system surfaces. Immersion-type disinfectants are not approved for use on system surfaces.
- DO NOT scratch or press on the panel with any sharp objects, such as pencils or pens, as this may result in damage to the LCD panel.
- Make sure not to spill or spray any liquid on the controls, into the system cabinet, or in the probe connection receptacle.
- Prior to cleaning, turn OFF power to the system and disconnect the mains cable.

Cleaning and Disinfecting the system



WARNING

All cleaners and disinfectants **NOT** on this list in *Table 3-8 Compatible Chemicals for Cleaning and disinfection* on page 199 are **unapproved** by GE HealthCare. Failure to follow guidelines could result in damage to the device.



CAUTION

When processing the operator control panel, make sure not to spill or spray any liquid on the controls, into the system cabinet, or in the probe connection receptacle.



CAUTION

To avoid liquids entering the product, **DO NOT** spray any liquid directly onto the surfaces. **ALWAYS** use a cloth or wipe.

These cleaners/disinfectants can be used on the console (Operator Panel, Monitor, Probe Holders, etc.), except for the probes. *Probe Intermediate-Level Disinfection - Wipe* on page 310 for probe disinfectant information and web links.

Always consult the cleaner or disinfectant manufacturer's instructions for proper use of their product. Wear appropriate PPE as indicated by the manufacturer.

Appropriate cleaners/disinfectants for the console that have been validated for compatibility are shown below:



WARNING

User shall read the instruction and be trained before doing cleaning and disinfection. Ensure that you follow the cleaning procedure and use the cleaning agents provided in this manual.

Table 3-8 Compatible Chemicals for Cleaning and disinfection

Chemical	System Cabinet	Operator Control Panel	Monitor and Touch Panel	Track Ball	Probe Holder and Gel Holder
PDI Sani-Cloth Plus	X	X	X	X	X
Protex	X	X	X	X	X
Tristel Wipes	X	X	X	X	X
Wet Wipe	X	X	X	X	X
Sono Ultrasound Wipes	X	X	X	X	X
Clorox wipes	X	X	X	X	X
Cavi Wipes	X	X	X	X	X
PDI Super Sani-Cloth Plus	X	X	X	X	X
PDI Easy Screen Cleaning Wipe	X	X	X	X	X
70% isopropyl alcohol	X	X	X	X	X
NOTE Effective Disinfection is always a balance between safe inactivation of infectious agents and undesirable side effects. Due to the generally uneven and irregular surface of Ultrasound consoles, a comprehensive surface disinfection process cannot be recommended by the manufacturer.					

Console

The ultrasound system console includes the System enclosure, Monitor, Operator Panel, Probe Holders and Gel Warmer. For Probes Reprocessing, see *Probe Intermediate-Level Disinfection - Wipe* on page 310.

Cleaning the Console

Always clean visible soil from surfaces first before disinfecting the Console.

Follow the cleaning/disinfecting frequency suggested in *Maintenance Schedule* on page 196.

To clean the system,

1. Moisten a soft, non-linting cloth with a mild, general purpose, non-abrasive soap and water solution or approved cleaning/disinfecting agent.

NOTE

The cloth/wipe should be damp, not dripping wet and running. Moisture should not drip into the crevices anywhere on the console.

NOTE

Refer to *Table 3-8 Compatible Chemicals for Cleaning and disinfection* on page 199 for a list of compatible solutions to be used on the Console.

2. Use a gentle wiping action to clean any surface on the console.

**HINTS**

A scrubbing action with the wipe may be necessary to help remove stubborn soil from the surfaces. However be careful with this action over cervix and gaps in the surface to prevent liquid from being scraped off the wipe and entering the product.

3. Wipe off excess cleaning agents.

NOTE

Do not spray any liquid directly into the unit.

NOTE

DO NOT scratch or press on the panel with any sharp objects, such as pencils or pens, as this may result in damage to the LCD panel.

Disinfecting the Console

For disinfectants to be effective, the surface must first be clean.

ALWAYS follow the manufacturer's instructions concerning the use of the disinfectant and follow the contact time to be sure the disinfectant's claims are accomplished.

Follow the cleaning/disinfecting frequency suggested in *Maintenance Schedule* on page 196.

Disinfect the desired surfaces of the console. To prevent cross contamination, surfaces that are often touched during exams should be disinfected after every patient. To disinfect the system,

1. Moisten a non-linting cloth with a liquid disinfectant or remove pre-moistened disinfectant wipe from the container.

**CAUTION**

If a cleaner/disinfectant wipe was used to clean off visible soil per the above section, a second, fresh cleaner/disinfectant product should be used for the disinfectant step.

2. Wet the surfaces by gently applying the cloth or wipe. Avoid high pressure or squeezing the wipe to avoid having the liquid enter the gaps and cervices of the Console. Scrubbing is not necessary in the disinfecting step; evenly applying the liquid is the goal.
3. Let the surface remain wet for the appropriate contact time.
4. If the surface does not remain wet for the full contact time, apply an additional application of the disinfectant, as necessary, to extend the time.
5. After the contact time has expired, remove excess liquid with a dry sterile cloth.
6. To avoid disinfectant buildup, or to remove disinfectant residue which may cause skin irritation, perform a rinse step with a sterile damp cloth.

Monitor

To clean the monitor (LCD Panel and monitor frame):

1. Moisten a soft, non-linting cloth with a mild, general purpose, non-abrasive soap and water solution.

NOTE

The cloth should be damp, not dripping wet.

2. Wipe down the top, front, back, and both sides of the monitor.
3. Wipe off excess cleaning agents with a non-linting cloth.

NOTE

Never use thinner, benzene, alcohol (ethanol, methanol, or isopropyl alcohol), abrasive cleaners, or other strong solvents, as these may cause damage to the monitor.

NOTE

DO NOT scratch or press on the LCD panel with any sharp objects, such as pencils or pens, as this may result in damage to the panel.

For disinfection for this part, refer to system level disinfection in *Disinfecting the Console* on page 200.

Control panel and keyboard

NOTE

Diligent cleaning of the console reduces the risk of spreading infection from person to person, and also helps to maintain a clean working environment.

Only use the following cleaners on the Control panel:

1. A non-abrasive soap and water solution (e.g. Palmolive Dishwashing Liquid, manufactured by Colgate-Palmolive).
2. Sani Wipes Alcohol-free (manufactured by Microgen Inc).
3. T-Spray II (manufactured by Pharmaceutical Innovations, Inc).

To clean the Control panel:

1. Turn off the power of the system.
2. Moisten a soft, non-linting cloth with water or a mild, non-abrasive soap and water solution.
3. Gently wipe the surface of the console.
4. Use a cotton swab to clean around keys or controls. Use a toothpick to remove solids between keys and controls.

When cleaning the operator control panel, make sure not to spill or spray any liquid on the controls, into the system cabinet, or in the probe connection receptacle.

5. In the event that disinfection is required or any stubborn stains remain, absorb a small quantity of isopropyl rubbing alcohol on a soft, dust-free cloth. Wipe the surface of the console. Make sure no liquid drips on or between the keys. Allow to dry.

To clean the keyboard:

1. Clean the keyboard as described above for the Control panel.

After The Exam Is Over

2. In the event that disinfection is required or any stubborn stains remain, absorb a small quantity of isopropyl rubbing alcohol on a soft, dust-free cloth. Wipe the surface of the key caps. Make sure no liquid drips on or between the keys. Allow to dry.

NOTE

When cleaning the operator control panel, make sure not to spill or spray any liquid on the controls, into the system cabinet, or in the probe connection receptacle.

For disinfection for this part, refer to system level disinfection in *Disinfecting the Console* on page 200.

Probe Cleaning Notes

Cleaning and disinfection solution for probe and console are different. Please avoid damaging the control panel while cleaning/disinfecting probes. Refer to the Probes Chapter, for probe cleaning and disinfecting instructions.



CAUTION

NEVER use any cleaner or disinfectant containing alcohol.

Probe shall not stay in probe holder on the Ultrasound system during cleaning/disinfection. If you use a spray cleaner, spray AWAY from the Ultrasound system. Spray can damage the controls.

Figure 3-22 DO NOT Spray a Probe While in its Holder



Cleaning the Trackball


Dust is often building up behind the ball, so it interferes with the ball rotation and for optical trackballs the light used for sensing. To get access for cleaning, you need to remove the ball.

The ball is held in position by the Dust Gasket.

1. Power off the system.

2. Rotate the retainer counterclockwise until it can be removed from the keyboard.

Figure 3-23 Rotate the retainer



3. Separate the trackball and the retainer. Wipe off any oil or dust from the trackball, retainer and the trackball housing using a cleaner or cotton swab.
4. Assemble the trackball and retainer, then put it into the housing and rotate it clockwise until its notches are set in position.



CAUTION

When cleaning, make sure not to spill or spray any liquid into the trackball housing (keyboard or system). DO NOT use cleaning agents containing ammonia, such as Windex, Screen-Clean, which will remove the anti-glare coating on the Trackball.

Probe holder

1. Remove the holder.

Figure 3-24 Remove the probe holder



2. Wash the holder with mild soap in lukewarm water. Scrub it using a soft non-linting sponge, gauze, or cloth to remove all visible residue from the surface. Prolonged soaking or scrubbing with a soft bristle brush (such as a toothbrush) may be necessary if material has dried onto the surface.

After The Exam Is Over

3. Rinse the holder with enough potable water.
4. Dry with a soft cloth and put it back.

For disinfection for this part, refer to system level disinfection in *Disinfecting the Console* on page 200.

System Cabinet

To clean the system cabinet:

1. Moisten a soft, non-linting cloth with a mild, general purpose, non-abrasive soap and water solution or a general purpose disinfectant.
2. Wipe down the top, front, back and both sides of the cabinet. Do not spray any liquid directly onto the unit.
3. In the event that disinfection is required or any stubborn stains remain, absorb a small quantity of isopropyl rubbing alcohol on a soft, dust-free cloth. Wipe the system cabinet and allow to dry.

NOTE

Do not spray any liquid directly into the unit.

Other Maintenance

Footswitch

To clean the footswitch:

1. Moisten a soft, non-abrasive folded cloth with a mild, general purpose, non-abrasive soap and water solution.
2. Wipe the external surfaces of the unit then dry with a soft, clean, cloth.

Cleaning the air filters

Clean the system's air filters to ensure that a clogged filter does not cause the system to overheat and reduce system performance and reliability. It is recommended the filters be cleaned every two weeks, but the requirements will vary due to your system use.



CAUTION

Be sure to lock the wheels before cleaning the air filters to avoid injury by any unexpected movement of the system.

DO NOT operate the unit without the air filters in place.

Allow the air filters to dry thoroughly before re-installing them on the unit.

1. Pull the air filter from the back cover of the system.

Figure 3-25 Remove the air filter



2. Dust the filters with a vacuum cleaner and/or wash it with a mild soapy solution.

NOTE

If washed, RINSE and DRY the filters before re-installation.

3. Put back the air filters.

Prevention of static electricity interference

Interference from static electricity can damage electronic components in the system. The following measures help to reduce the likelihood of electrostatic discharge:

- Wipe the alphanumeric keyboard and monitor with lint-free tissue or a soft cloth dampened with anti-static spray on a monthly basis.
- Spray carpets with anti-static spray because constant walking on carpets in or near the scanning room may be a source of static electricity.

Disposal

Table 3-9 WEEE symbol

Symbol	Location
	Rear of the system and Probe connector





This symbol indicates that waste electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately. Please contact an authorized

representative of the manufacturer for information concerning the decommissioning of your equipment.

Troubleshooting

Please refer to Service Manual if other messages appear on the monitor display.

Table 3-10 Error messages

 Warning The system has detected the lower air filter requires cleaning. Please clean the lower air filter. OK	The system has detected the lower air filter requires cleaning. Please clean the lower air filter. 1. Shutdown the system. Disconnected the power cord. 2. Clean the air filter according to <i>Cleaning the air filters</i> on page 204.
 Warning System temperature is too high. System will shut down. Current Progress: 6% Ok Cancel	System temperature is too high. System will shut down. 1. Shutdown the system. 2. Clean the air filter according to <i>Cleaning the air filters</i> on page 204.
 Warning System voltage fault. System will shut down. Current Progress: 6% Ok Cancel	System voltage fault. System will shut down. 1. Select OK and reboot the system. 2. If the same message appears after reboot, shut down the system and turn off the breaker. Then turn on the system according to <i>Power On</i> on page 44.
 Error System Error. Please Restart the system. Ok	System Error. Please reboot the system. 1. Select OK and reboot the system 2. If the same message appears after reboot, shut down the system and turn off the breaker. Then turn on the system according to <i>Power On</i> on page 44.


InSite ExC

InSite ExC is your direct link with a GE Online Service Engineer or Applications Support Engineer or a Request for Service via the InSite ExC  link at the upper right of the monitor.

Types of InSite ExC Service

1. **ContactGE.** Opens a service dispatch with GE HealthCare Service.
2. **Start/Stop Agent.** Enable/Disable Agent.
3. **ConnectToGE.** Agent is set to fast poll after running the command.
4. **Connect Clinical Lifeline.** Agent is set to fast poll, turn on disruptive mode, if disruptive mode is off. Turn on VNC, if VNC server is not running.
5. **EZ Configuration Wizard.** Use to configure some common system settings.
6. **Temperature And Power.** Use to check the temperature log.
7. **EZ Diagnostic.** Use to analyze the system error automatically.
8. **Versana Club.** Use to enter Versana Club.
9. **Disruptive and VCO.** Use to enable Disruptive mode and open VCO. Click it and the following window will appear.

Figure 3-26 Select Enable



Here is the explanation for current design on "Disruptive mode and VCO" button function:

- a. The button is designed to Enable VCO in quick way, so when clicking "Enable", it will enable both "Disruptive mode" and "VCO" (Disruptive mode ON is the prerequisites for VCO ON).
 - b. When clicking "Disable", it can only disable "VCO", but will not disable the Disruptive mode.
 - c. The user has to enter CSD to disable the Disruptive mode if he/she wants to do that.
10. **EZ Settings Restore.** Use to quick backup and restore network configuration, user defined settings and RSVP configuration.
 11. **Download.** Use to download the latest software package.
 12. **Cancel.** Use to abort the downloading task and clear the downloaded data.

Initiating a Request for Service (RFS)

To initiate an RFS,

After The Exam Is Over

1. Position the Windows pointer on top of the GE InSite ExC icon at the top right of the display.
2. Press **Set** Key, select **ContactGE**. This opens of the RFS screen which sends a service dispatch directly to GE HealthCare Service after you fill in the following information:
 - Problem Type
 - Problem Area
 - Problem Description
 - First Name
 - Last Name
 - Phone Number
 - Email
3. After you have completed filling in all of this information, press **Submit** to initiate the Request for Service.

After the request is submit successfully, the following status appears:

Figure 3-27 Request for Service is submitted



All requests for service are listed on the **RFS History** for your review. You can review your RFS details by clicking **Details** on RFS History list.

Figure 3-28 RFS History

Problem Type	Problem Area	First Name	Last Name	Status	Date	Link
Service	Hardware Electrical	test	test	Success	2019-05-06 14:32:27 -07:00	Details ^
Application	Software	Li	Yuanyuan	Success	2019-05-21 10:04:02 -07:00	Details
Service	Software	Tingting	Chen	Failed	2019-05-22 13:59:13 -07:00	Details
Service	Software	gsd	sefgs	Success	2019-05-22 14:06:09 -07:00	Details
Service	Software	Tina	Chen	Success	2019-05-22 14:09:29 -07:00	Details v

Initiating a Technical or Clinical Support Request

To initiate Technical or Clinical Support,

1. Position the Windows pointer on top of the GE HealthCare InSite ExC icon at the upper right of the display.
2. Press the **Set** Key.

NOTE

When you have contacted the Online Center/Field Engineer (OLC/FE), you may be asked to click on **ContactGE** to increase the polling rate so that the OLC/FE can connect more quickly.

3. Select **ConnectToGE** for Technical Support or press **Cancel** to exit.

NOTE

In addition to contacting a technical/clinical support person, selecting this also changes the polling time from 15 minutes to 15 seconds so that your call can be answered as quickly as possible, as well as allowing disruptive mode.

InSite ExC icons appear differently, depending on their state:

Table 3-11 InSite Icons

	RSvP Agent is not registered and Disruptive mode is off.
	RSvP Agent is not registered and Disruptive mode is on.
	RSvP Agent is connected and Disruptive mode is off.
	RSvP Agent is connected and Disruptive mode is on.
	Device is remotely connected.

InSite ExC Definitions

Here are definitions for the different InSite ExC states:

Virtual Console Observation (VCO). Allows Technical Support to control the ultrasound system functionality remotely.

Disruptive. Allows GE HealthCare's Technical Support person to connect to your system via VCO, to run diagnostics directly on your ultrasound system, and to collect system logs. When the system is in Disruptive Mode, the icons are red. There are two disruptive states. If you see a telephone with a clock, then the system is in Disruptive, Not Connected Mode. If you see a telephone with GE HealthCare, then the system is in Disruptive, Connected Mode.

Non-Disruptive. Allows GE's Technical support person to look around on your system, but cannot perform any service-related functions, depending on whether InSite has connected or not connected. There are two Non-Disruptive states. If you see a black and white icon, InSite ExC is activated, but not open for Technical Support access. If you see a yellow icon, InSite ExC is activated and the Technical Support person can look around on your system, but cannot perform any service-related functions.

Connected. InSite ExC is connected.

Not Connected. InSite ExC is not connected

NOTE

When Disruptive mode has been activated or a diagnostic has been run, the message, "Due to Service testing reboot required," appears in red at the bottom of the display. It is recommended that you reboot the system before use. Make sure you disable disruptive mode before rebooting or the message will not be cleared.

Exiting InSite ExC

To exit InSite ExC,

1. Exit InSite ExC page.
2. Reboot your ultrasound system.

EZ Configuration Wizard


EZ Configuration wizard is a function to enable the operator to configure some common system settings when turning on the system for the first time after the software installation.

NOTE

Password setting is required when turning on the system for the first time after the software installation.

User can also enter EZ Configuration Wizard by clicking the Insite icon at the top right of the screen.

Figure 3-29 Enter EZ Configuration Wizard



Then follow below steps to complete configuration.

1. Select password level in **Password Policies** or **Skip For Now** to setup later. User can change the policy level later by **Utility > Admin > User policies**. Users can select the option **Lowest**, **Medium**, **High** or **Highest** for policy level. The complexity of the password will be differed by the policy level you selected.


NOTE

Password policy level setting is only required when turning on the system for the first time or after the software installation.

NOTE

When user select lowest, a warning will display to inform that the current security policy will put your system at risk.

Figure 3-30 Select password policy level

**NOTE**

If user select lowest level in password policies, ADM password can be empty during Admin Registration.

NOTE

Exit button is not activated in password policies setting, users are not able to exit the EZ Configuration Wizard in this step.

- Set ADM password in **Admin Registration**.

NOTE

Admin Registration is only required when turning on the system for the first time or after the software installation.

NOTE

Exit button is not activated in Admin Registration setting, users are not able to exit the EZ Configuration Wizard in this step.

NOTE

ADM password can be empty if password policy level is selected as LOWEST.

NOTE

A valid password must be at least 8 characters long and has a maximum length of 256 characters.

Figure 3-31 ADM password setting

Password Policies

Admin Registration

Language Setting

General Information

Network Information

Summary Information

Create ADM Account

Operator: ADM

Enter Password:

Confirm password:

Auto Logon: **Enabled**

Minimum Password Requirements:

- No password complexity rules

The passwords are critical for your system. Please store the passwords safely.

Previous Next Exit

If users did not set ADM password during system setup. User can also set operator password later by **Utility > Admin > Users**.

3. Press **Next** to select the appropriate language for system language and keyboard language from the drop-down list.

Figure 3-32 System Language settings

Password Policies

Admin Registration

Language Setting

General Information

Network Information

Summary Information

Language Setting

Language: ENG

Keyboard: ENG_US

NOTE: If you want to see the following pages in the selected language, please restart the system.

Previous Next Exit


NOTE

All required setting after turning on system for the first time is done. If you select **Exit**, you will be able to exit the EZ Configuration Wizard and start using your system.

- If you do not change the language, press Next to continue.
- If you change the language setting, a language setup window will display.

Set language and keyboard language, then press **OK** to restart the system.

Figure 3-33 Language settings confirmation



NOTE

If you press **Previous**, you will go to the previous page.

4. This screen shows the hospital and time information, you can set the system date and time on this page.

Figure 3-34 General Information

The image shows a "General Information" screen with a sidebar on the left containing a menu with the following items: Password Policies, Admin Registration, Language Setting, General Information (highlighted), Network Information, and Summary Information. The main content area is divided into two sections: "Hospital Information" and "Time Information".

Hospital Information

- Hospital: Hospital Name
- Department: Development
- Customer Name:
- Address:
- Continent: AFRICA
- Country: ALGERIA

Time Information

- Time Zone: (UTC-08:00) Pacific Time (US & Canada)
- System Date: A calendar for August 2024 is displayed. The date 30 is selected.
- System Time: 07:21:45. A checkbox for "24 Hour" is checked.

At the bottom of the screen are three buttons: "Previous", "Next", and "Exit".

Press **Next** to continue.

After The Exam Is Over

5. The **Network Information** screen shows the configuration of wireless and local network.

Figure 3-35 Network Information

The screenshot shows the 'Network Information' configuration screen. On the left is a sidebar menu with options: Password Policies, Admin Registration, Language Setting, General Information, Network Information (highlighted), and Summary Information. The main area is titled 'Choose The Connection' with a dropdown menu set to 'Ethernet(Connected)'. Below this is the 'TCP/IP' section with a checked 'Enable DHCP' option. The configuration fields are as follows:

Field	Value
IP Address	10 . 191 . 63 . 102
Subnet Mask	255 . 255 . 255 . 0
Default Gateway	10 . 191 . 63 . 1
DNS	10 . 56 . 56 . 56

At the bottom of the TCP/IP section is an 'Apply' button. At the bottom right of the screen are three buttons: 'Previous', 'Next', and 'Exit'.

Press **Next** to continue.

6. **Summary Information** shows the report of the previous settings. User can export it to the database.

Figure 3-36 Summary Information


The screenshot shows the 'Summary Information' screen. The sidebar menu is the same as in Figure 3-35, with 'Summary Information' highlighted. The main area displays a summary of the settings under the heading 'IWSummary'. The summary is organized into sections:

- general_info**
 - hospital : Hospital Name
 - department : Development
 - customer :
 - address :
 - timezone : (UTC-08:00) Pacific Time (US & Canada)
 - continent : AFRICA
 - country : ALGERIA
- lang_setting**
 - lang : ENG
 - keyboard :
- network_setting**

At the bottom right is an 'Export' button. At the bottom of the screen are three buttons: 'Previous', 'Next', and 'Finish'.

Press **Export** and click **OK** to save the summary information.

Figure 3-37 Export Summary



7. Press **Finish** to exit **EZ Configuration Wizard**. The system will display a warning to restart to activate your configuration. Press **OK** to restart.

NOTE

System will restart only after user configures **EZ Configuration Wizard** for the first time or new software version is installed.

NOTE

For better security, it is recommended to change the passwords of these two accounts when using the device for the first time. Go to **Utility > Service > Utilities > Change Password**, change the passwords of **GEService** and **SVCservice**. The complexity of the new password cannot be lower than the password policy level that has been set (**Lowest**, **Medium**, **High** or **Highest**).

eDelivery - Software update

External software environment of the ultrasound system should be updated regularly for better use experience. If necessary, the external software environment of the ultrasound system can also be updated urgently.

User can update to the latest software in two ways.

As part of the product lifecycle management, GE HealthCare regularly analyzes and integrates software updates from our third party vendors into our products. These are typically released as part of regular updates or software releases.

The two available downloading options for eDelivery are:

- Through the GE HealthCare service platform on the ultrasound system. This requires Insite RSvP connectivity.
- Through an end-user portal to a local storage location (i.e. a readable / writeable flash drive with enough storage space) and install it on the ultrasound system.

Software update via Insite Remote Service Platform (RSvP)

Software update for the system may become available for download and installation through the GE HealthCare Service platform.

Users must have administrator rights to perform the software download and installation. A user who is not logged in as ADM (administrator) will see the notification of an available update, but not be allowed to initiate the download.



CAUTION

Please backup presets and patient data when updates system software. Remote software download should not change user presets or affect customer database; however, it is always best practice to ensure patient data and preset are backed up before proceeding with any software installation.

NOTE

Be sure to end any open exams before software update.

NOTE

Please allow approximately half an hour for complete software download (the download time may vary due to network connection speed).

NOTE

Please allow approximately 30 minutes for complete installation.

NOTE

Software upgrade through the GE HealthCare service platform may not be available in all markets.

Software download and installation




CAUTION

Please end current exam before downloading the package provided by eDelivery.

1. Press the notification icon to display the service menu, click **Download** to check the update.


Figure 3-38 Service Notification



2. Available software updates are displayed in the list. If you want to refresh the query for available updates, press **Refresh**.


3. Select the desired software version and press **Download** to download the software package, or press **Exit** to exit the window.

Figure 3-39 Download or Exit



4. During the software download process, the status will be displayed as "**Downloading**". You can press **Pause** to suspend the download or press **Cancel** to abort the download.

Figure 3-40 Downloading

**NOTE**


While the software is downloading, the **Exit** button is disabled. If you want to go back to scanning page, please press **Pause** to suspend the download process first and then press the **Exit** button.

- When the download is in progress, press the **Pause** button to suspend the download process.
- When the download process is suspended, press the **Resume** button to recover the download process from the point where it is stopped.
- If there is an error during the downloading process, press the **Retry** button to recover the download process from the beginning.

NOTE

Before installation, user should login with the ADM account, otherwise there will be pop-up message to remind. Only after logging in with the ADM account, the user will have access to install software.

Figure 3-41 Access Denied



- When the download process finishes, the software is ready to be installed in the ultrasound system. Log in with the ADM account, then click **Yes** to start the installation.

Figure 3-42 Install to the system





Figure 3-43 Start Installation



NOTE

Please allow enough time for ultrasound device to reboot automatically to complete the installation process.

- The system will reboot several times to complete the installation.

NOTE


DO NOT power off the system during the software installation.

NOTE

A typical full installation may take up to 30 minutes.

7. After software installation is complete and the system is rebooted, a **Reminder after Installation** window will appear. This is to remind the user to manually restore configurations if the system has related settings backed up before the software update.

Figure 3-44 Reminder after Installation



8. Click **OK** to confirm the installation is completed.

NOTE

Ensure that the entire system functions normally after service procedures.

Software update via End-User Portal

Customers who have purchased the end-user eDelivery option will receive a customer account to download software within the End-User Portal.

Users are created for the account based on e-mail addresses provided by the customer at the point of sale. These e-mail addresses are the log-in credential for the End-User Portal along with a temporary password provided to the user through e-mail. When logging in to the End-User Portal the first time, the user is prompted to change the password and enter a secret question and answer for password retrieval.

Follow the below instructions to download software from the portal:

1. Log on to the portal website which is provided to end user via a welcome email:
<https://gehealthcare.flexnetoperations.com/flexnet/operationsportal>

2. Log in using the user name (e-mail) and password.

Figure 3-45 Login Window



The login window features the GE HealthCare logo at the top left. Below the logo is a "Login" section with three input fields: "Username", "Password", and a language dropdown menu currently set to "English (United States)". A "Forgot password?" link is located below the password field. A "Log in" button is positioned to the right of the language dropdown.

3. The **License & Delivery Portal** dashboard is displayed. All downloadable packages will be listed in a link under **Your Downloads**.


Figure 3-46 Your Downloads



After The Exam Is Over


- Click the link to enter the **Download Packages** page. Downloadable packages with links will be listed. Click on the desired link to enter the **Downloads** page.

Figure 3-47 Select link



- Right click on the link under **File Name**, select **Save link as...** to download the file to your storage device (a readable and writeable USB flash drive with enough capacity).

Figure 3-48 Save Link as...



NOTE

Multiple files may be downloaded on one device.

- Load the downloaded software on the ultrasound system from the selected storage location.

Loading the Software

To load a software onto the ultrasound system,


1. Power off the ultrasound system.
2. Insert the USB Flash Drive into a USB port.

NOTE

Ensure that the system is USB Device Enabled.

3. Power on the ultrasound system. The software program files will be recognized automatically, select **Install SW...** on the **Start Application** screen.

Figure 3-49 Select Install SW...

**NOTE**

If by accident you try to load a software that is not compatible with the software on the ultrasound system, an error message will indicate "The package present in media is not compatible. Please contact GE HealthCare Service" and "Software installation is not started".

4. Click **OK** on **StartLoader** screen. Then select the package you want to install, and click **INSTALL** to start installation.
5. When the installation is complete, the system will reboot automatically.
6. Remove the USB stick from the ultrasound device.
7. The whole installation is completed once system reboots.

NOTE

Installation time depends on the type of update. No action is needed until the installation is completed.

Reload the Software

Reload Software provides specific software versions for users to reload.

NOTE


Users are required to be logged into the system as administrators (ADM) in order to be granted access.

To reload the software onto the ultrasound system:

After The Exam Is Over


1. Press **Utility** > **System** > **About**, click **Reload Software** to enter the operation page.

Figure 3-50 Reload Software



2. Select the software version needed, then click **OK**.

Figure 3-51 Select the software version




- By default, current installed software version will be displayed, user can reload current software version to rollback software changes to the initial installed state.
- After a new full software package is installed, current version will be displayed as previous version, user can select to reload current version or previous version.

NOTE

After rollback to previous software version, the "Current/Previous" reload options will alternate.

3. A pop-up message will appear to remind user to continue or cancel the reload process. Click **Yes** to start reloading software or **No** to cancel.

Figure 3-52 Pop-up message



Assistance

Supplies/Accessories



CAUTION

DO NOT connect any probes or accessories without approval by GE HealthCare.



CAUTION

Use only GE HealthCare approved internal equipment when replacing an internal peripheral. The user or the operator should never install/replace the internal peripheral. Service representatives authorized by GE HealthCare will install/replace the internal peripheral.

Not all features, products, probes or peripherals described in this document may be available or cleared for sale in all markets. Please contact your local GE HealthCare Ultrasound representative to get the latest information.

Contact the distributor, GE HealthCare affiliate or sales representative for approved peripherals. For HCATs, contact your sales person.

The following supplies/accessories have been verified to be compatible with the system:

Peripherals

Table 3-12 Peripherals and Accessories

Accessory	Units
Printer USB isolation	Each
Gel Warmer	Each
Video Output Adapter	Each
HDMI cable option	Each
Wireless Adaptor	Each
Bluetooth Adaptor	Each
Transcend TS8XDVDS-K DVDRW kit	Each
UP-D898DC Printer	Each
BW PRINTER (UP-D898MD)	Each
SONY UPD25 Color Printer	Each
898 printer paper	Each
HP Officejet 200 Mobile series Printer	Each
USB Stick	Each
1TB mobile USB HDD	Each
Footswitch MKF 2 1S/1S-MED USB GP26	Each
1-Peadl type footswitch 'Whanam FSU-1000'	Each
USB barcode reader	Each
ECG ASSY	Each
1 TB SSD	Each
Digital Expert Connect	Each
Rechargeable Li-ion Battery Pack	Each
Power Cord	Each
Push-Push A/N keyboard	Each
ADULT TEE CLIP-ON BITE GUARD	Each
ADULT TEE CLIP-ON BITE GUARD OPR.	Each
ADULT TEE SCANHEAD PROTECTION COVER	Each
ADULT TEE CONVENTIONAL BITE GUARD	Each
BITE HOLE INDICATOR	Each
TEE STORAGE RACK	Each
China Control panel protect film	Each
A/N Keyboard Film	Each

NOTE

Other commercial printers not listed in above table must connect with PIT (Printer Isolation Transformer) when being used.

NOTE

The transformer does not support the printer with power more than 150VA.

NOTE

The ultrasound system supports two types printer isolation transformers: 110V/220V. Please use the appropriate type according to the local power supply.

Probe**Table 3-13 Probes and Accessories**

Probe	HCAT	Biopsy Guide	Biopsy Guide HCAT	Biopsy Guide Image
4C-RS Convex	H4000SR	Multi-angle, disposable needle guide with a reusable plastic bracket	E8385NA	
8C-RS Convex	H40402LS	Not Available	N/A	N/A
E8C-RS Convex	H40402LN	Fixed angle, disposable needle guide	E8385MJ	
		Fixed angle, reusable with a stainless steel bracket	H40412LN	
E8Cs-RS Convex	H48062AF	Fixed angle, disposable needle guide	E8385MJ	
		Fixed angle, reusable with a stainless steel bracket	H40412LN	
C1-5-RS Convex	H40462AL	Multi-angle, disposable needle guide with a reusable plastic bracket	H40432LE	

Probe	HCAT	Biopsy Guide	Biopsy Guide HCAT	Biopsy Guide Image
IC9-RS Convex	H48691PJ	Fixed angle, disposable needle guide	H48691YX	
		Fixed angle, disposable needle guide (with cover)	H48691YW	
		Fixed angle, reusable with a stainless steel bracket	H48701MN	
3Sc-RS Sector	H45041DL	Multi-angle, disposable needle guide with a reusable plastic bracket	H46222LC	
6S-RS Sector	H45021RP	Not Available	N/A	N/A
12S-RS Sector	H44901AB	Not Available	N/A	N/A
L6-12-RS Linear	H48062AC	Multi-angle (in plane biopsy kit), disposable needle guide with a reusable plastic bracket	H40432LC	
12L-RS Linear	H40402LY	Multi-angle (in plane biopsy kit), disposable needle guide with a reusable plastic bracket	H40432LC	
		Infinite-angle (in plane biopsy kit), disposable needle guide with a reusable plastic bracket	H48392LT	
		Transverse bracket (out of plane biopsy kit), disposable needle guide with a reusable plastic bracket	H48392LL	

Probe	HCAT	Biopsy Guide	Biopsy Guide HCAT	Biopsy Guide Image
L3-12-RS Linear	H44901AP	Multi-angle Bracket	H48302AA	
LK760-RS Linear	H44901AF	Not Available	N/A	N/A
L8-18i-RS Linear	H40462LF	Not Available	N/A	N/A
9L-RS Linear	H40442LL	Multi-angle, disposable needle guide, with a reusable plastic bracket	H4906BK	
L4-20t-RS Linear	H48062AJ	Verza biopsy starter kit	H45201BLF	
RAB2-6-RS 4D	H48681WR	Multi-angle, disposable needle guide with a reusable plastic bracket	H48681ML	
RIC5-9A-RS 4D	H48701EJ	Fixed-angle, disposable needle guide with a disposable plastic bracket	H48691Z	
		Fixed-angle, disposable needle guide with a disposable plastic bracket	H48681GF	
		Fixed-angle, reusable needle guide with a stainless steel bracket	H46721R	

Probe	HCAT	Biopsy Guide	Biopsy Guide HCAT	Biopsy Guide Image
E7C8L-RS Bi-plane	H48062AL	Fixed-angle, disposable needle guide with a reusable plastic bracket	H40202E	
		Fixed-angle, disposable needle guide	NA	NA
BE9CS-RS Bi-plane	H40482LN	Fixed angle, disposable needle guide	E8387M/ E8013AW	
		Fixed angle, reusable needle guide with a stainless steel bracket	E8387MA/ E8013AX	
6Tc-RS TEE	H45551ZE	Not Available	NA	NA

System Data

Features/Specifications

Table 3-14 Physical Attributes

Dimensions and Weight - Height: less than 1840 mm - Width: less than 600 mm - Depth: less than 970 mm - Weight: less than 73 kg without any probes or peripherals Console Design 3, 4 or 5 active probe ports - SSD (Solid State Disk) - Integrated speakers - Cable hooks - Handle - Removable air filter - Lever for moving control panel up/down and swivel - Gel Holder and Probe holder - Fixed or Flexible arm - Locking mechanism that provides caster lock	Keyboard - Alphanumeric keyboard - Interactive back-lighting - Digital keyboard Monitor - High-Resolution monitor - Brightness and contrast adjustment - LCD Fold-down and lock mechanism for transportation Electrical Power - Voltage: 100-240 V AC ($\pm 10\%$) - Frequency: 50/60 Hz - Power: Consumption less than 500VA with on board peripherals
--	--

Table 3-15 Inputs and Outputs Signal

External Inputs and Outputs	Outputs
• USB port • Ethernet	• Composite Video out (via an adaptor) • S-Video (via an adaptor) • HDMI video

Table 3-16 System Overview - Operating Modes and Peripheral Options

Operating Modes	Peripheral Options
• B-Mode • CHI-Mode • M-Mode • Color Flow Mode (CF) • PW (Pulse Wave) Doppler (with High PRF) • CW (Continued Wave) Doppler (option) • Combined Mode: B+B, B+M, B+CF, B+PDI, B+PW, B+CW, B+CF+PW, B+CF+CW, B+PDI+PW, B+PDI+CW, CHI+CHI, CHI+M, CHI+CF, CHI+PDI, CHI+PW, CHI+CW, CHI+CF+PW, CHI+CF+CW, CHI+PDI+PW, CHI+PDI+CW Note: Whether the combined modes listed above can be implemented depends on whether the corresponding independent mode options are installed. • Power Doppler Imaging (PDI) • Directional Power Doppler Imaging (Directional PDI) • Color M Mode (option) • Anatomical M-Mode (AMM) (option) • Curved Anatomical M-Mode (Curved AMM) (option) • Tissue Velocity Imaging (TVI) (option) • TVM (Tissue Velocity M Mode image) (option) • Tissue Velocity Doppler (TVD) (option) • Static 3D/Realtime 4D (option) • B Flow (option) • B Flow Color (option)	• Printer USB isolation • Gel Warmer • Video Output Adapter • HDMI cable option • Wireless adapter • Bluetooth adaptor • DVDRW kit • B&W (black and white) image printer • Color image printer • HP Officejet 200 Mobile series Printer • USB memory stick • 1TB mobile USB HDD • Footswitch • Barcode reader • ECG kit • Digital Expert Connect • Rechargeable Li-ion Battery Pack • Power Cord • Push-Push A/N keyboard • Report printer • External USB printer • Removable storage media • Isolation USB port • A/N keyboard film • SSA (Secured Service Access) • Second SSD • Endo-probe holder

Table 3-17 System Overview - System Features

System Features

• Contrast Image (option) • Elastography(option) • Easy 3D (option) • Advanced 3D (option) • TUI (Tomographic Ultrasound Imaging) combined in 4D (option) • Tissue Velocity Imaging (TVI) with QAnalysis (option) (option) • LOGIQView (option) • CrossXBeam (option) • SRI HD • Auto-IMT (option) • Stress Echo (option) • Auto EF (Ejection Fraction) (option) • Virtual Convex • Sonobiometry • VOCAL (option) • Breast productivity (including BI-RADS) (option) • Thyroid productivity (including TI-RADS) (option) • Breast Productivity-AI (option) • Thyroid productivity-AI (option) • Auto Bladder Volume measurement • Needle recognition (option) • Auto TGC • Auto Tissue Optimization (ATO) • Automatic Spectrum Optimization (ASO) • Whizz • Whizz CF	• Whizz Easy Style (option) • Digital LGC • Digital TGC • Digital Keyboard • V-Live 2.0 (option) • Strain Imaging (option) • Strain Rate (option) • Whizz Follicle (option) • Whizz RenderLive (option) • Breast Care (option) • Follow-up tool (option) • LI-RADS(ACR) (option) • ECG (option) • Cine • DICOM (option) • Scan Coach • Scan Assistant • My Trainer • Report/Worksheet • Multi-touch • Fast boot-Standby • Wide dual screen measurement • Synchronous mode • Voice Comment • VCO for customer (option) • Probe Check • Imaging Insights (option) • Whizz Report • Whizz Note
---	---

Quality Assurance

A good Quality Assurance Evaluation program consists of periodic systematic actions that provide the user with adequate confidence that their diagnostic ultrasound system will produce consistently high quality images and quantitative information.

Therefore, it is in the best interests of every ultrasound user to routinely monitor equipment performance.

The frequency of Quality Assurance evaluations should be based on user's specific needs and clinical practice.

Periodic monitoring is essential in order to detect the performance changes that occur through normal aging of system components. Routine equipment evaluations may also reduce the duration of exams, number of repeat exams, and maintenance time required.

For details on system and peripheral routine preventive maintenance instructions, see **System Care and Maintenance** section for more information.

Typical Tests to Perform

Quality assurance measurements provide results relating to system performance. Typically these are:

- Axial Measurement Accuracy
- Lateral Measurement Accuracy
- Axial and Lateral Resolution
- Penetration
- Functional & Contrast Resolution
- Gray Scale Photography.

With these tests, a performance baseline can be set at installation with the phantom in your department. Future test results can be compared to the baseline in order to maintain a record of system performance trends.

The phantom shown is shown as a representative example of a phantom. You can select from any number of phantoms available on the market.

Frequency of tests

Quality assurance tests are used to determine whether a scanner is providing the same level of performance from day to day.

The frequency of testing varies with the amount of system usage and modes to be tested. It is recommended that the user perform quality assurance tests at least every three months or every 400 patient studies. Tests should also be performed when a question about system performance exists.

A mobile system may require more frequent tests.

Image quality should also be tested immediately after the following events:

- Service calls
- System upgrades/modifications
- Dropped probe, power surge, etc.


Phantoms

Quality Assurance Evaluations may be done with phantoms and test objects that are applicable to the parameters being evaluated or to the user's clinical practice.

Typical phantoms are composed of material that acoustically mimic human tissue. Pins, anechoic and echogenic targets are physically positioned to provide information for a variety of tests.

The RMI 403GS phantom is shown in the illustration below as a representative example of a phantom.

Figure 3-53 Representative Phantom Example



1. Penetration
2. Axial Distance Measurement
3. Functional Resolution
4. Lateral Resolution
5. Lateral Distance Measurement
6. Axial Resolution
7. Contrast Resolution and Gray Scale Photography
8. Gray Scale Plane Targets

Baselines

An absolute necessity for a quality assurance program is establishing baselines for each test or check. Baselines are established after the system has been verified to be working properly at installation or after a repair. If a probe or major assembly is replaced, new baselines should be generated.

Baselines can be made by adjusting system parameters to prescribed levels or to the best possible image. The key factor to remember is reproducibility. The same conditions must be reproduced for each periodic check.

All system parameters not displayed on the monitor should be recorded for the permanent record.

Periodic Checks

Periodic checks should be performed in accordance with your facility's quality assurance requirements. For the data to be valid, periodic checks should mimic the baseline setup parameters.

The resulting image, when scanning the phantom exactly as before, should be recorded and compared to the baseline. When a matching image is obtained, it can be assumed that the system performance has not degraded from the baseline.

If a significant difference between the baseline and periodic check is noted, double check the system setup and repeat the test. If the difference between the baseline and periodic check persists, contact a local Service Representative.

Failing to reproduce the control settings as in the baselines will introduce errors in the data and potentially invalidate the results.

Results

Lack of standardization among test instruments, the wide range of acceptance criteria, and incomplete knowledge regarding the significance of certain performance parameters prohibit the establishment of absolute performance criteria for these tests.

Quality Assurance Evaluation results should be compared to previously-recorded results.

Performance trends can then be detected. Unacceptable performance or diminishing trends should be identified for maintenance or repair before a malfunction or inappropriate diagnosis occurs.

The user should determine the best method for recording and archiving the baseline and periodic checks. In most cases the choice is hard copy.

It is important to maintain good consistent records for inspections that may arise, as well as to detect system performance trends.

System Setup

The user should tailor the tests to their particular needs. It is certainly not necessary to make all checks with all probes. A representative example, with the probes used most often by the customer, should be adequate in judging system performance trends.

Use a gray scale phantom as the scan object for the tests. Commercial phantoms are supplied with its own operator manual. Be familiar with proper phantom operating procedures prior to use for quality assurance evaluations.

1. Adjust image monitor. Brightness and contrast should be set to the normal viewing of a good gray scale image.
2. Check all recording devices for proper duplication of image monitor. Ensure that what is seen is what is recorded.

After The Exam Is Over

3. Annotate non-displayed image processing controls.
4. Set TGC slide pots to center (detent) position.
5. Place focal zone marker(s) in area of interest for an optimum image.

Test Procedures

The following are recommended Quality Assurance tests. A brief description of the test, the benefit it provides and steps to accomplish the test are supplied.

The importance of recording scan parameters and consistent record keeping cannot be stressed enough. Reproducibility to monitor system trends is the key to quality assurance evaluations.

Using the system's dual image display format is often very convenient and saves recording media.

Axial distance measurements

Description

Axial measurements are the distance measurements obtained along the sound beam. See *Figure 3-53 Representative Phantom Example* on page 234 for more information.

Benefit

The accurate measurement of the size, depth and volume of a structure is a critical factor in determining a proper diagnosis. Most imaging systems use depth markers and/or electronic calipers for this purpose.

Method

Axial distance should be measured in the near, mid and far fields as well as in zoom. If necessary, different depths or fields of view can be tested.

Procedure

To measure axial distance:

1. Scan a test phantom with precisely-spaced vertical pin targets. Adjust all scan controls, as necessary, for the best image of the pin targets to typical depths for the probe being used.
2. Press **Freeze** to stop image acquisition and perform a standard distance measurement between the pins at different points in the image. Record all images for archiving.
3. Scan the vertical pins in zoom or at different depth/scale factors.
4. Press **Freeze** to stop image acquisition; repeat the distance measurements between pins and record the images for archiving.
5. Document the measurements for reference and future comparison.

Contact a Service Engineer if vertical measurements differ by more than 1.50% of the actual distance.

Lateral distance measurements

Description

Lateral measurements are distance measurements obtained perpendicular to the axis of the sound beam.

Benefit

The purpose is the same as vertical measurements. Precisely-spaced horizontal pin targets are scanned and results compared to the known distance in the phantom.

Method

Lateral distance should be measured in the near, mid and far fields as well as in zoom. If necessary, different depths of fields of view can be tested.

Procedure

To measure lateral distance:

1. Scan a test phantom with precisely-spaced horizontal pin targets. Adjust all scan controls, as necessary, for the best image of the pin targets from side to side.
2. Press **Freeze** to stop image acquisition and perform a standard distance measurement between the pins at different points in the image. Record all images for archiving.
3. Scan the horizontal pins in zoom or at different depth/scale factors.
4. Press **Freeze** to stop image acquisition; repeat the distance measurements between pins and record the images for archiving.
5. Document the measurements for reference and future comparison.

Contact a Service Engineer if horizontal measurements differ by more than 3mm or 3% of that distance, whichever is greater.

Axial resolution

Description

Axial resolution is the minimum reflector separation between two closely-spaced objects to produce discrete reflections along the axis of the sound beam. It can also be monitored by checking the vertical size of known pin targets.

Benefit

In clinical imaging, poor axial resolution displays small structures lying close together as a single dot. This may lead to improper interpretation of the ultrasound image.

Procedure

To measure Axial resolution:

1. Scan a test phantom with precisely-spaced vertical pin targets.
2. Adjust all scan controls, as necessary, for the best image of the pin targets to typical depths for the probe being used.
3. Press **Freeze** to stop image acquisition.

After The Exam Is Over

4. Perform a standard distance measurement of the pin vertical thickness at different points in the image. Record all images for archiving.
5. Scan the vertical pins in zoom or at different depth/scale factors.
6. Press **Freeze** to stop image acquisition; repeat the vertical thickness measurements of the pins and record the images for archiving.
7. Document the measurements for reference and future comparison.

Axial resolution should remain stable over time. Contact a Service Engineer if any changes are observed.

Lateral resolution

Description

Lateral resolution is the minimum reflector separation between two closely spaced objects to produce discrete reflections perpendicular to the axis of the sound beam. It can also be monitored by checking the horizontal size of known pin targets. See *Figure 3-53 Representative Phantom Example* on page 234 for more information.

Lateral resolution is dependent upon the beam width produced by the probe. The narrower the beam, the better the lateral resolution.

The beam width is affected by the frequency, degree of focusing, and distance of the object from the face of the probe.

Benefit

Clinically, poor lateral resolution will display small structures lying close together as a single dot. This may lead to improper interpretation of the ultrasound image.

Procedure

To measure lateral resolution:

1. Scan a test phantom with precisely-spaced horizontal pin targets.
2. Adjust all scan controls, as necessary, for the best image of the pin targets from side to side.
3. Press **Freeze** to stop image acquisition and perform a standard distance measurement of the horizontal thickness of a pin at different points in the image. Record all images for archiving.
4. Scan the horizontal pins in zoom or at different depth/scale factors.
5. Press **Freeze** to stop image acquisition; repeat the horizontal thickness measurements of the pins and record the images for archiving.
6. Document the measurements for reference and future comparison.

Pin width should remain relatively constant over time ("1mm). Dramatic changes in pin width may indicate beamforming problems. Contact a Service Engineer if beam width changes consistently over 2 to 3 periodic tests.

Penetration

Description

Penetration is the ability of an imaging system to detect and display weak echoes from small objects at large depths. See *Phantoms* on page 234 for more information.

Penetration can be affected by the system's:

- Transmitter/receiver
- Degree of probe focusing
- Attenuation of the medium
- Depth and shape of reflecting object
- Electromagnetic interference from local surroundings.

Benefit

Weak reflecting echoes are commonly produced from the internal structure of organs. Definition of this tissue texture is important in the interpretation of the ultrasound findings.

Method

Scan a phantom to see how echoes begin to fade as depth is increased. The maximum depth of penetration is the point at which homogeneous material in the phantom begins to lose brightness.

Procedure

To measure penetration:

1. Set the front panel TGC slide pots to their center (detent) position.
2. Gain and acoustic output can be adjusted, as necessary, since these values are displayed on the monitor.
3. Scan a test phantom along the vertical pin targets to typical depths for the probe being used.
4. Perform a standard distance measurement from the top of the image displayed to the point at which homogeneous material in the phantom begins to lose brightness.
5. Document the depth measurement for reference and future comparison.

Contact a Service Engineer if the depth of penetration shifts more than one centimeter (1cm) when using the same probe and same system settings.

Functional resolution

Description

Functional resolution is an imaging system's ability to detect and display the size, shape, and depth of an anechoic structure, as opposed to a pin target.

The very best possible image is somewhat less important than reproducibility and stability over time. Routine tests at the same settings should produce the same results.

Benefit

The data obtained will give a relative indication of the smallest structure the system is capable of resolving at a given depth.

Procedure

To measure functional resolution:

1. Set the front panel TGC slide pots to their center (detent) position.

After The Exam Is Over

2. Gain and acoustic output can be adjusted as necessary, since these values are displayed on the monitor.
3. Scan a test phantom with a vertical row of anechoic cyst targets to typical depths for the probe being used.
4. Evaluate the cysts at various depths for a good (round) shape, well-defined borders and no fill in. Remember, TGC slide pots are centered and should remain fixed. This may NOT provide optimal cystic clearing.
5. Document all results for future reference and comparison.

Contact a Service Engineer if a greatly distorted image is obtained.

Contrast resolution

Description

Contrast resolution is the ability of an imaging system to detect and display the shape and echogenic characteristics of a structure.

Specific values measured are less important than stability over time. Routine tests at the same settings should produce the same results.

Benefit

A correct diagnosis is dependent upon an imaging system's ability to differentiate between a cystic or solid structure versus echo patterns from normal surrounding tissue.

Method

A phantom with echogenic targets of different sizes and depths should be used.

Procedure

To measure contrast resolution:

1. Set the front panel TGC slide pots to their center (detent) position. Set dynamic range to 54 db.
2. Gain and acoustic output can be adjusted, as necessary, since these values are displayed on the monitor.
3. Scan a test phantom with echogenic targets at the depths available.
4. Evaluate the echogenic targets for contrast between each other and between the surrounding phantom material. Remember, TGC slide pots are centered and should remain fixed. This may NOT provide an optimal scan image.
5. Document all results for future reference and comparison.

Contact a Service Engineer if the echogenic characteristics or shapes of the targets appear distorted.

Gray Scale photography

Description

Poor photography will cause loss of low level echoes and the lack of contrast between large amplitude echoes.

Benefit

When photographic controls and film processors are properly adjusted, weak echoes, as well as strong echoes, are accurately recorded on film.

Procedure

1. Adjust the camera according to the manufacturer's instructions until the hard copy and video display are equal.
2. Scan the phantom and it's echogenic contrast targets.
3. Make a hard copy photograph of the display and compare it to the image on the video monitor for contrast and weak echo display.
4. Document all results for future reference and comparison.

Contact a Service Engineer if camera cannot duplicate what is on the image monitor.

NOTE

Optimization of brightness/contrast controls on the display monitor is imperative in order to make sure that the hardcopy and monitor look alike.

The display monitor is adjusted first. The hardcopy camera or printer is adjusted to match the display monitor.

Setting up a Record Keeping System

Preparation

The following is needed:

- Quality Assurance binder.
- Hard copy or electronic file of images.
- Quality Assurance Checklists.
- Display the following information while testing quality assurance:
 - Acoustic Output
 - Gain
 - Depth
 - Probe
 - Dynamic Range
 - Set up new patient to be the name of the test.
- Annotate the following:
 - Any control where its value is NOT displayed.
 - Significant phantom information.

Record Keeping

Complete the following:

After The Exam Is Over

1. Fill out the Ultrasound Quality Assurance Checklist for each probe, as scheduled.
2. Make a hard copy or archive the image.
3. Compare images to baseline images and acceptable values.
4. Evaluate trends over previous test periods.
5. File hard copy or electronic file of images and checklist in Quality Assurance binder.

Ultrasound Quality Assurance Checklist

Table 3-18 Ultrasound Quality Assurance Checklist (Part 1)

Performed By		Date
System		Serial Number
Probe Type	Probe Model	Serial Number
Phantom Model	Serial Number	Room Temperature
Acoustic Output	Gain	Focal Zone
Gray Map	TGC	Depth
Monitor Setting		
Peripheral Settings		
Other Image Processing Control Settings		

Table 3-19 (Part 2)

Test	Baseline Value Range	Tested Value	Image Hardcopy/ Archived	Acceptable ? Yes/No	Service Called (Date)	Date Resolved
Vertical Measurement Accuracy						
Horizontal Measurement Accuracy						
Axial Resolution						
Lateral Resolution						
Penetration						
Functional Resolution						
Contrast Resolution						
Gray Scale Photography						

Probe Check

Probe Check is a probe assessment tool that evaluates each probe element. Probes have to be clean, any gel residue will provide incorrect data. This test is an assessment that is intended to be used comparatively during the life of the probe to evaluate possible probe

deterioration over time. The system default setting is running Probe Check every time for USA region with software option control.



CAUTION

Probe Check can detect most probe defects but not ALL possible probe defects.



CAUTION

GE HealthCare is NOT responsible for the confirmation about the normal function of the ultrasound probe before use.

Probe check will evaluate the probe performance and provide the user with an indication of potential impact to the diagnostic image if it is compromised due to transducer malfunction.



CAUTION

DO NOT allow the probe head to hang free. Impact to the probe head could result in irreparable damage.



CAUTION

DO NOT unplug the probe during Probe Check process, otherwise the test result shall be failed.

Pass/Fail criteria

NOTE

If any one of the two Pass or Fail criteria listed below are met, the result is a failure.

1. Total number of element failures exceeded.
2. Number of adjacent element failures exceeded (based on a pre-defined gap between elements which are considered adjacent).

Supported probes

Probe Check supports on all probes for this system.

Presets

Probe Check provides five interval options for auto triggering, press **Utility > System > General > Probe Check** to select the preset.

- **Never:** probe checking will never automatically opened.
- **Every time:** probe checking will be automatically opened each time while triggering probe active.
- **Once 1 day:** probe checking will be automatically opened one time every day while triggering probe active.
- **Once 7 days:** probe checking will be automatically opened one time every week while triggering probe active.
- **Once 30 days:** probe checking will be automatically opened one time every month while triggering probe active.

NOTE

The system default is every time for USA region. For other regions, the system default is Never. The first boot up shall follow this setting.




CAUTION

User shall perform probe check manually if the auto trigger interval is not configured as "every time" when changing a probe with the same name as previous probe.

NOTE

For all options other than NEVER, Probe Check will be triggered automatically while switching the probes. Switching the presets individually will not trigger Probe Check.

Figure 3-54 Probe Check preset



Probe Test Auto Trigger

NOTE

Always remember to clean the probe surface before probe check procedure.

1. Check that the probe to be tested is thoroughly clean and dry. For cleaning the probe, refer to *Care and Maintenance* on page 300 for more information. Connect it to the probe port on the scanner, then hold the probe in the air ready for testing.

NOTE

The auto trigger assessment is only available when the **Probe Check** preset is **Every Time**.

NOTE

The auto trigger assessment is only available for the probe connected to the activated probe port.

NOTE


Prior to inserting the probe, ensure that the connector locking handle is positioned to the unlock state.

NOTE

Ensure the probe is locked to the system before you start the probe diagnostics.


2. Press **Probe** on touch panel, select the probe then the system will do probe test automatically, a message will pop-up on the monitor "**Probe Check: Testing probe, please wait...**".
3. If the probe test is Passed, there will be a message "**Probe Check: Probe test passed.**" showing on the monitor.

Figure 3-55 Probe Check: Probe test passed



4. If the probe test is Failed, the checking results will be shown on the pop-up window.

Figure 3-56 Probe Test Fail



NOTE

The testing results may vary in different series of GE HealthCare ultrasound products due to differentiated design. If you have any concern with the testing result, please contact with your service partner.

5. Follow the suggestions shown on the pop-up window to contact with your service partner, follow the guidance from your service partner to determine if there is issue with the probe or the system itself.

Manual Trigger


A manual probe check should be conducted when images acquired with a deficient probe look abnormal. Phenomena include, but are not limited to overall decreased brightness or/and partial signal loss.

After The Exam Is Over

Example phantom scanning images are provided below for your reference. Whenever you see similar cases or are not confident in the status of a probe, please perform probe check manually.


- Condition 1: Total number of element failures exceeded.

Figure 3-57 Abnormal image acquired with a deficient L6-12-RS probe that has distributed defect elements



- Condition 2: Number of adjacent element failures exceeded.


Figure 3-58 Abnormal image acquired with a deficient L6-12-RS probe that has adjacent defect elements



Users can perform probe check manually:

1. Press **Probe** on touch panel, select the probe to be tested, click **Test Probe** button.

Figure 3-59 Manual Probe Test on Touch Panel



2. Click **Test Probe** button on the pop-up window. The testing workflow is same as step 2 to step 5 in **Auto Trigger**.

Figure 3-60 Manual Probe Test




Image Quality Check (IQC)

Image Quality Check is intended to facilitate Image Quality checks during Quality Assurance Evaluations. Quality Assurance tests are used to determine whether a scanner is providing the same level of performance year after year.

By using the same settings year after year, this ensures that the data collection is consistent, independently of who performs the test.

This preset only includes fundamental settings for B-Mode. Processing modes like SRI, Harmonics, etc., are turned off.

NOTE

A tissue-mimicking phantom should be used for IQC evaluations.

To performn an Image Quality Check (IQC):

1. Activate IQC via **Utility > Imaging Preset Manager > Category** (select the Category first).
2. Click on the plus sign in front of IQC for Service and select **IQC**.
3. Assign IQC to a Touch Panel key by using the right arrow key.
4. Map the IQC to the location you want it to appear on the Touch Panel.
5. Select **Probe**. Then select **IQC**.

Privacy and Security

Potential Hazardous Situations Resulting from Failures of the IT Network

The following general hazardous situations have been identified as potentially hazardous as a result of the IT network failing to provide the required characteristics specified above.

- Delayed or impaired access to images or other exam information or patient data.
- Permanent loss of images or other exam information or patient data.
- Corruption of images or other exam information or patient data.
- LAN connection loss during operation may cause loss of data, and damage data integrity.

Warning

In addition to the hazardous situations identified above, connection of the Ultrasound System Software to a network that includes other equipment could result in other unidentified risks to patients, operators or third parties. The responsible organization should identify, analyze, evaluate and control these risks on an ongoing basis including after changes to the network such as those listed below, which could introduce new risks and require additional analysis.

- Changes in network configuration.
- Connection of additional items to the network.
- Disconnecting items from the network.
- Update of equipment connected to the network.
- Upgrade of equipment connected to the network.

Contact Information

Contacting GE HealthCare Ultrasound

For additional information or assistance, please contact your local distributor or the appropriate support resource listed on the following pages:

Internet

<http://www.gehealthcare.com>

<https://cleaning.gehealthcare.com>

Clinical Questions

For information in the United States, Canada, Mexico and parts of the Caribbean, call the Customer Answer Center. TEL: (1) 800-682-5327 or (1) 262-524-5698

In other locations, contact your local Applications, Sales, or Service Representative.

Service Questions

For service in the United States, call GE HealthCare CARES.

TEL: (1) 800-437-1171

In other locations, contact your local Service Representative.

Information Requests

To request technical product information in the United States, call GE HealthCare.

TEL: (1) 800-643-6439

In other locations, contact your local Applications, Sales, or Service Representative.

Placing an Order

To order accessories, supplies, or service parts in the United States, call the GE HealthCare Technologies Contact Center.

TEL: (1) 800-558-5102

In other locations, contact your local Applications, Sales, or Service Representative.

Table 3-20 Americas

ARGENTINA	GE HealthCare Argentina Nicolas de Vedia 3616 piso 5 Buenos Aires - 1307	TEL: (+54) 11-5298-2200
BRAZIL	GE HealthCare do Brasil Comércio e Serviços para Equipamentos Médicos - Hospitalares Ltda. Av. Magalhães de Castro, 4800, Andar 10 Conj. 101 e 102, Torre 3 Cidade Jardim CEP: 05676-120 - São Paulo/SP - Brasil	TEL: 3004 2525 (Capitals and Metropolitan Regions) 08000 165 799 (Other Locations)
CANADA	GE HealthCare Ultrasound 9900 Innovation Drive Wauwatosa, WI 53226	TEL: (1) 800-668-0732 Customer Answer Center TEL: (1) 262-524-5698

LATIN & SOUTH AMERICA	GE HealthCare Ultrasound 9900 Innovation Drive Wauwatosa, WI 53226	TEL: (1) 262-524-5300 Customer Answer Center TEL: (1) 262-524-5698
MEXICO	GE HealthCare Sistemas Medicos de Mexico S.A. de C.V. Rio Lerma #302, 1° y 2° Pisos Colonia Cauhtemoc 06500-Mexico, D.F.	TEL: (5) 228-9600 FAX: (5) 211-4631
USA	GE HealthCare Ultrasound 9900 Innovation Drive Wauwatosa, WI 53226	TEL: (1) 800-437-1171 FAX: (1) 414-721-3865

Table 3-21 Asia

ASIA PACIFIC JAPAN	GE HealthCare Asia Pacific 4-7-127, Asahigaoka Hinoshi, Tokyo 191-8503, Japan	TEL: +81 42 585 5111
AUSTRALIA	32 Phillip Street Parramatta 2150 Sydney, NSW, Australia	TEL: 1800 659 465
CHINA	GE HealthCare - Asia No. 1, Yongchang North Road Beijing Economic & Technology Development Area Beijing 100176, China	TEL: (8610) 5806 8888 FAX: (8610) 6787 1162 Service: 4008128188 (24h)
INDIA	Wipro GE HealthCare Pvt Ltd No. 4, Kadugodi Industrial Area Sadaramangala, Whitefield Bangalore, 560067	TEL: +(91) 1-800-425-8025
KOREA	15F, 416 Hangang Dae ro, Chung-gu Seoul 04637, Korea	TEL: +82 2 6201 3114
NEW ZEALAND	Level 7 Vero Centre 48 Shortland St, Auckland, 1010 New Zealand	TEL: 0800 659 465
SINGAPORE	GE HealthCare ASEAN (Singapore) 11 North Buona Vista Drive #11-07 The Metropolis Tower 2 Singapore 138589	TEL: +65 6291 8528 FAX: +65 6291 7006

Table 3-22 Africa

EGYPT	GE Medical Systems Egypt, LLC Plot 44 Tesseen El Shamaly Street Al Salam Axis First Sector City Centre, 5th settlement Cairo, Egypt	TEL: +20 2 25354200 FAX: +20 2 25370031
KENYA	GE HealthCare East Africa Services Limited General Mathenge Drive, Courtyard Building Westlands Nairobi 30 00100 KE	TEL: +254 719 093 044
NIGERIA	GE HealthCare International Operations (Nig) Ltd Bishop Aboyade Cole Street No. 927/928 Mansard Place, PO Box 54255 Victoria Island Lagos LA NG	TEL: +234 (01) 4607101 TEL: +234 (01) 4607102
KINGDOM OF SAUDI ARABIA	GE HealthCare Arabia Co. Ltd Platinum Centre, Building 1 Salahuddin Ayoubi Road Riyadh-12811 Kingdom of Saudi Arabia	TEL: +966 (11) 494 5779 FAX: +966 (11) 207 3946
SOUTH AFRICA	GE HealthCare South Africa (Pty) Ltd. 60 Glenhove Road Green on Glenhove Customer Innovation Centre Johannesburg GP 2196 ZA	TEL: +270100725000 FAX: +27 0862958385

Table 3-23 Europe

AUSTRIA	GE HealthCare Austria GmbH & Co OG EURO PLAZA, Gebäude E Technologiestrasse 10 A-1120 Vienna	TEL: (+43) 1 97272 0 FAX: (+43) 1 97272 2222
BELGIUM & LUXEMBURG	GE HealthCare BVBA/SPRL Kouterveldstraat 20 1831 DIEGEM	TEL: (+32) 2 719 7204 FAX: (+32) 2 719 7205
CZECH REPUBLIC	GE Medical Systems Česká Republika, s.r.o. Bucharova 2641/14 158 00 Praha 5 Česká republika	TEL: (+420) 224 446 162 FAX: (+420) 224 446 161
DENMARK	GE HealthCare Park Allè 295 DK-2605 Brøndby, Denmark	TEL: (+45) 43 295 400 FAX: (+45) 43 295 399

ESTONIA & FINLAND	GE HealthCare Finland Oy Kuortaneenkatu 2, 000510 Helsinki P.O. Box 330, 00031 GE Finland	TEL: (+358) 10 39 48 220 FAX: (+358) 10 39 48 221
FRANCE	GE Medical Systems SCS Division Ultrasound 24 Avenue de l'Europe - CS20529 78457 Vélizy Villacoublay Cedex	TEL: (+33) 1 34 49 52 70 FAX: (+33) 13 44 95 202
GERMANY	GE HealthCare GmbH Beethovenstrasse 239 42655 Solingen	TEL: (+49) 212-28 02-0 FAX: (+49) 212-28 02-380
GREECE	GE HealthCare 8-10 Sorou Str. Marousi Athens 15125 Hellas	TEL: (+30) 210 8930600 FAX: (+30) 210 9625931
HUNGARY	GE HealthCare Hungary Zft. Bence utca 3 Budapest BU 1138 HU	TEL: (+36) 23 410 314 FAX: (+36) 23 410 390
IRELAND	NORTHERN IRELAND GE HealthCare Victoria Business Park 9, Westbank Road Belfast BT3 9JL. REPUBLIC OF IRELAND GE HealthCare 3050 Lake Drive Citywest Business Campus Dublin 24	TEL: (+44) 028 90229900 TEL: 1800 460 550 FAX: (+353) 1 686 5327
ITALY	GE Medical Systems Italia spa Via Galeno, 36, 20126 Milano	TEL: (+39) 02 2600 1111 FAX: (+39) 02 2600 1417
KAZAKHSTAN	«Дженерал Электрик Қазақстан» ЖШС Қазақстан, Алматы қаласы, Медеу ауданы, көшесі ЗЕНКОВ, үй 26/41, пошталық индексі 050010	T +7 727 3560020
LUXEMBORG	See Belgium.	
NETHERLANDS	GE HealthCare De Wel 18 B, 3871 MV Hoevelaken PO Box 22, 3870 CA Hoevelaken	TEL: (+31) 33 254 1290 FAX: (+31) 33 254 1292
NORWAY	GE HealthCare Vingmed Ultrasound AS Sandakerveien 100C 0484 Oslo, Norway GE HealthCare Vingmed Ultrasound Strandpromenaden 45 P.O. Box 141, 3191 Horten	TEL: (+47) 23 18 50 50 FAX: (+47) 23 18 60 35 TEL: (+47) 33 02 11 16

After The Exam Is Over

POLAND	GE Medical Systems Polska Sp. z o.o., ul. Woloska 9 02-583 Warszawa, Poland	TEL: (+48) 22 330 83 00 FAX: (+48) 22 330 83 83
PORTUGAL	GE HealthCare Portuguesa SA Avenida do Forte 6 - 6A Edificio Ramazzotti 2790-072 Carnaxide	TEL: (+351) 21 425 1300 FAX: (+351) 21 425 1343
RUSSIA	GE HealthCare Presnenskaya nab. 10 Block C, 12 floor 123317 Moscow, Russia	TEL: (+7) 4957 396931 FAX: (+7) 4957 396932
SPAIN	GE HealthCare España C/ Gobelos 35-37 28023 Madrid	TEL: (+34) 91 663 2500 FAX: (+34) 91 663 2501
SWEDEN	GE HealthCare Sverige AB FE 314, 182 82 Stockholm Besöksadr: Venevagen 89 Danderyd, Sverige	TEL: (+46) 08 559 500 10 FAX: (+46) 08 559 500 15 Service Center (+46) 020-120 14 36
SWITZERLAND	GE Medical Systems (Schweiz) AG Europastrasse 31 8152 Glattbrugg	TEL: (+41) 1 809 92 92 FAX: (+41) 1 809 92 22
TURKEY	GE HealthCare Türkiye Istanbul Office Levent Ofis Esentepe Mah. Harman Sok. No:8 Sisli-Istanbul	TEL: +90 212 398 07 00 FAX: +90 212 284 67 00
UNITED ARAB EMIRATES (UAE)	GE HealthCare Dubai Internet City, Building No. 18 First Floor, Dubai - UAE	TEL: (+971) 4 429 6101 or 4 429 6161 FAX: (+971) 4 429 6201
UNITED KINGDOM	GE Medical Systems Ltd Pollards Wood Nightingales Lane Chalfont St Giles Buckinghamshire HP8 4SP	TEL: (+44) 1494 544000 FAX: (+44) 1707 289742
For all other European countries not listed, please contact your local GE HealthCare distributor or the appropriate support resource listed on www.gehealthcare.com .		

Manufacturer



GE Medical Systems (China) Co., Ltd.
No. 19, Changjiang Road
Wuxi National Hi-Tech Development Zone
214028 Jiangsu
P.R. China

TEL: +86 510 85225888; FAX: +86 510 85226688

Factory Sites

GE Medical Systems (China) Co., Ltd.
No. 19, Changjiang Road
Wuxi National Hi-Tech Development Zone
214028 Jiangsu
P.R. China

The product from below factory site can only be available in India:

Wipro GE Medical Device Manufacturing Private Limited
No. 4 Kadugodi Industrial Area, Sadarmangala, Whitefield
Bangalore, Karnataka, 560067 INDIA

Chapter 4

Safety

In this section

Owner Responsibility.	258
Bioeffect Disclosure.	259
Safety Precautions.	260
Device Labels.	287

Owner Responsibility

Owner requirements

It is the responsibility of the owner to ensure that anyone operating the system reads and understands this section of the manual. However, there is no representation that the act of reading this manual renders the reader qualified to operate, inspect, test, align, calibrate, troubleshoot, repair or modify the system. The owner should make certain that only properly trained, fully-qualified service personnel undertake the installation, maintenance, troubleshooting, calibration and repair of the equipment.

The owner of the ultrasound unit should ensure that only properly trained, fully-qualified personnel are authorized to operate the system. Before authorizing anyone to operate the system, it should be verified that the person has read, and fully understands, the operating instructions contained in this manual. It is advisable to maintain a list of authorized operators.

Should the system fail to operate correctly, or if the unit does not respond to the commands described in this manual, the operator should contact the nearest field GE HealthCare Ultrasound Service Office.

For information about specific requirements and regulations applicable to the use of electronic medical equipment, consult the local, state and federal agencies.

Notice against user modification

Never modify this product, including system components, software, cables, and so on. User modification may cause safety hazards and degradation in system performance. All modification must be done by a GE HealthCare qualified person.

Bioeffect Disclosure

Following are potential risks inherent to technology:

1. During a diagnostic ultrasound examination, high frequency sound penetrates and interacts with tissue in and around the area of anatomy to be imaged. Only a small portion of this sound energy is reflected back to the transducer for use in constructing the image while the remainder is dissipated within the tissue. The interaction of sound energy with tissue at sufficiently high levels can produce biological effects (aka bioeffects) of either a mechanical or thermal nature. Although the generation of bioeffect is intentional with therapeutic ultrasound, it is generally undesired in diagnostic application and may be harmful in some conditions.
2. Cavitation may occur by Ultrasound technology during a diagnostic ultrasound examination. Mechanical Index (MI) is an attempt to indicate the probability that cavitation might occur within the tissue. Always observe the Acoustic Output display for possible effects.

Training

During each ultrasound examination the user is expected to weigh the medical benefit of the diagnostic information that would be obtained against the risk of potentially harmful effects. Once an optimal image is achieved, the need for increasing acoustic output or prolonging the exposure cannot be justified. It is recommended that all users receive proper training in applications before performing them in a clinical setting. Contact the GE HealthCare sales representative for training assistance.

ALARA



CAUTION

- Ultrasound can produce harmful effects in tissue and potentially result in patient injury.
- Always minimize exposure time and keep ultrasound levels low when there is no medical benefit. Use the principle of ALARA (As Low As Reasonably Achievable), increasing output only when needed to obtain diagnostic image quality. Observe the acoustic output display and be familiar with all controls affecting the output level. See the Bioeffects section of the Acoustic Output chapter in the Advanced Reference Manual for more information.



CAUTION

- Ultrasound can produce harmful effects in tissue and potentially result in patient injury.
- The operator of the device must sufficiently understand the acoustic output and be able to obtain the related thermal index values. Always minimize exposure time to the irradiation and keep ultrasound acoustic output level low for embryos or fetuses.

Safety Precautions

Precaution Levels

Various levels of safety precautions can be found on the equipment. Different levels of concern are identified by one of the following flag words and icons, which precede the precautionary statement.



DANGER

A danger symbol and icon indicates that a specific hazard is known to exist which, through inappropriate conditions or actions, will cause :

- Severe or fatal personal injury



WARNING

A warning symbol and icon indicates that a specific hazard is known to exist which, through inappropriate conditions or actions, may cause:

- Minor personal injury
- Substantial property damage (substantial property damage is defined as damage requiring service for the system to function)



CAUTION

A caution symbol and icon indicates that a potential hazard may exist which, through inappropriate conditions or actions, will or can cause:

- Property damage

NOTE

A note indicates precautions or recommendations that should be used in the operation of the ultrasound system, specifically:

- Maintaining an optimum system environment
- Using this Manual
- Notes to emphasize or clarify a point

Hazard Symbols

Icon Description

Potential hazards are indicated by the following icons:

Table 4-1 Potential Hazards

Icon	Potential Hazard	Usage
	Biological Hazard Describes precautions necessary to prevent the risk of disease transmission or infections. • Patient/user infection due to contaminated equipment.	• Cleaning and care instructions • Sheath and glove guidelines
	Electrical Hazard Describes precautions necessary to prevent the risk of injury through electric hazards. • Electrical micro-shock to patient, e.g., ventricular.	• Probes • ECG, if applicable • Connections to back panel
	Moving Hazard Describes precautions necessary to prevent the risk of injury through moving or tipping hazard! • Console, accessories or optional storage devices that can fall on patient, user, or others. • Collision with persons or objects may result in injury while maneuvering or during system transport. • Injury to user from moving the console.	• Moving • Using brakes • Transporting
	Acoustic Output Hazard • Patient injury or tissue damage from ultrasound radiation.	• ALARA, the use of Power Output following the 'as low as reasonably achievable' principle
	Explosion Hazard Describes precautions necessary to prevent the risk of injury through explosion hazard. • Risk of explosion if used in the presence of flammable anesthetics.	• Flammable anesthetic
	Fire and Smoke Hazard • Patient/user injury or adverse reaction from fire or smoke. • Patient/user injury from explosion and fire.	• Replacing fuses • Outlet guidelines

Important Safety Considerations

The following topic headings (Patient Safety, and Equipment and Personnel Safety) are intended to make the equipment user aware of particular hazards associated with the use of this equipment and the extent to which injury can occur if precautions are not observed. Additional precautions may be provided throughout the manual.



DANGER

Improper use can result in serious injury. The use of the system outside the described conditions or intended use, and disregarding safety related information is considered abnormal use. The user must be thoroughly familiar with the instructions and potential hazards involving ultrasound examination before attempting to use the device. Training assistance is available from GE HealthCare if needed.

Disregarding information on safety is considered abnormal use.

Following are potential risks inherent to technology:

- Ultrasonic energy delivered to non-targeted tissue with the use of Ultrasound devices, and the interaction of sound energy with tissue at sufficiently high levels can produce biological effects.
- Monitoring the Mechanical Index (MI) can be a tool to help monitor the probability that cavitation could occur.

Reporting

Any serious incident related to the use of this GE HealthCare ultrasound device should be reported to both the Manufacturer and the health authority/competent authority where the device is installed.

To report to GE HealthCare:

Either contact your local service representative or report to:

in-box.complaints@gehealthcare.com.

Please provide the following information:

- The catalogue number or the model designation of the device as stated on its identification plate affixed on the device.
- The SystemID/serial number/lot number of the device.
- Date of incident.
- Description of incident, including any patient or user impact/injury.
- Your contact information(facility, address, contact name, title, and telephone number).

Patient Safety Related Hazards



WARNING

The concerns listed can affect the safety of patients undergoing a diagnostic ultrasound examination.

Patient Identification

Always include proper identification with all patient data and verify the accuracy of the patient's name and ID numbers when entering such data. Make sure correct patient ID is provided on all recorded data and hard copy prints. Identification errors could result in an incorrect diagnosis.

The ultrasound system is not meant to be long term storage for patient data or images. The customers are responsible for the data on the system and a regular backup is highly recommended.

It is advisable to back up system data prior to any service repairs to the hard drive. It is always possible during system failure and repair to lose patient data. GE HealthCare will not be held responsible for the loss of this data.

Diagnostic Information

The images and calculations provided by the system are intended for use by competent users, as a diagnostic tool. They are explicitly not to be regarded as the sole, irrefutable basis for clinical diagnosis. Users are encouraged to study the literature and reach their own professional conclusions regarding the clinical utility of the system.

The user should be aware of the product specifications and of the system accuracy and stability limitations. These limitations must be considered before making any decision based on quantitative values. If in doubt, the nearest GE HealthCare Ultrasound Service Office should be consulted.

Equipment malfunction or incorrect settings can result in measurement errors or failure to detect details within the image. The equipment user must become thoroughly familiar with the equipment operation in order to optimize its performance and recognize possible malfunctions. Applications training is available through the local GE HealthCare representative. Added confidence in the equipment operation can be gained by establishing a quality assurance program.



WARNING

The system provides calculations (e.g. estimated fetal weight) and charts based on published scientific literature. The selection of the appropriate chart and clinical interpretation of calculations and charts is the sole responsibility of the user. The authorized user should consider proper indications for the use of a calculation or chart as described in the scientific literature. The diagnosis, decision for further examination, and medical treatment must be performed by qualified personnel following good clinical practice.



CAUTION

Features that facilitate measurements such as VOCAL or SonoNT must be used with extreme care. The measurement results are a suggestion of the system, if in doubt verify with manual measurement methods.

The user is responsible for the diagnostic interpretation of the measurement results.

Mechanical hazards

The use of damaged probes or improper use and manipulation of intracavity probes can result in injury or increased risk of infection. Inspect probes often for sharp, pointed, or rough surface damage that could cause injury or tear protective barriers. Become familiar with all instructions and precautions provided with special purpose probes.



ELECTRICAL HAZARD

A damaged probe can also increase the risk of electric shock if conductive solutions come in contact with internal live parts. Inspect probes often for cracks or openings in the housing and holes in and around the acoustic lens or other damage that could allow liquid entry. Become familiar with the probe's use and care precautions outlined in Probes and Biopsy.



DANGER

Ultrasound transducers are sensitive instruments which can easily be damaged by rough handling. Take extra care not to drop transducers and avoid contact with sharp or abrasive surfaces. A damaged housing, lens or cable can result in patient injury or serious impairment or operation.



WARNING

Observe probe immersion levels.

Inspect probes for sharp edges or rough surfaces that could injure sensitive tissue.

DO NOT bend or pull the cable forcefully, to avoid mechanical shock or impact to the probe.

Scanner and electrosurgical units



CAUTION

This equipment provides no special means of protection from high frequency (HF) burns that may result from using an electrosurgical unit (ESU). To reduce the risk of HF burns, avoid contact between the patient and ultrasound transducer or ECG electrodes while operating the ESU. Where contact cannot be avoided, make sure the transducer or ECG electrodes are not located between the ESU active and dispersive electrodes and keep the ESU cables away from the transducer or ECG cables.

Equipment and Personnel Safety

The concerns listed in the Equipment and Personnel Safety section can seriously affect the safety of equipment and personnel during a diagnostic ultrasound examination.

Equipment and Personal Safety Related Hazards



DANGER

This equipment contains dangerous voltages that are capable of serious injury or death.

If any defects are observed or malfunctions occur, stop operating the equipment and perform the proper action for the patient. Inform a qualified service person and contact a Service Representative for information.

There are no user replaceable or serviceable components inside the console. Refer all servicing to qualified service personnel only.

Ensure that unauthorized personnel do not tamper with the unit.



ELECTRICAL HAZARD

To avoid injury:

- Do not remove protective covers. No user serviceable parts are inside. Refer servicing to qualified service personnel.
- To assure adequate grounding, connect the attachment plug to a reliable (hospital grade) grounding outlet (having equalization conductor ).
- Never use any adaptor or converter of a three-prong-to-two-prong type to connect with a mains power plug. The protective earth connection will loosen.
- Do not place liquids on or above the console. Spilled liquid may contact live parts and increase the risk of shock.
- In North America, a 220-240V installation requires the use of a center-tapped AC power source.

**WARNING**

Any electronic device can fail without warning signs, therefore the user is advised to follow local clinical practice guidelines for having a backup imaging plan when performing time-critical image-guided examinations and interventions.

**WARNING**

Remove probe from patient body before defibrillation application.

**WARNING**

Keep the system away from magnetic resonance imaging (MRI) equipment.

**WARNING**

The internal circuits of the unit use high voltages, capable of causing serious injury or death by electrical shock.

**WARNING**

To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

**CAUTION**

In case any serious incident occurred in relation to this ultrasound system should be reported to GE Healthcare and the competent Authority.

**DANGER**

The concerns listed below can seriously affect the safety of equipment and personnel during a diagnostic ultrasound examination.

**EXPLOSION HAZARD**

Risk of explosion if used in the presence of flammable anesthetics.

Never operate the equipment in the presence of flammable or explosive liquids, vapors or gases. Malfunctions in the unit, or sparks generated by fan motors, can electrically ignite these substances. Operators should be aware of the following points to prevent such explosion hazards.

- If flammable substances are detected in the environment, do not plug in or turn on the system.
- If flammable substances are detected after the system has been turned on, do not attempt to turn off the unit, or to unplug it.
- If flammable substances are detected, evacuate and ventilate the area before turning off the unit.

**SMOKE AND FIRE HAZARD**

The system must be supplied from an adequately rated electrical circuit. The capacity of the supply circuit must be as specified.

**DANGER****Biological Hazard**

For patient and personnel safety, be aware of biological hazards while performing invasive procedures. To avoid the risk of disease transmission:

- Use protective barriers (gloves and probe sheaths) whenever possible. Follow sterile procedures when appropriate.
- Thoroughly clean probes and reusable accessories after each patient examination and disinfect or sterilize as needed. Refer to Probes and Biopsy for probe use and care instructions.
- Follow all infection control policies established by your office, department or institution as they apply to personnel, equipment and accessories.

**DANGER**

Do not unpack the ultrasound system. This must be performed by qualified service personnel only. Improper unpacking could lead to injury.

**DANGER**

GE HealthCare recommends dedicated probes for use on humans only or animals only. Mark probes dedicated for animals with special labels.

Observe any country specific rules and regulations for handling equipment used on both animals and humans. Such national restrictions may prohibit transfer of probes used on animals to humans and vice-versa.

Failure to follow these instructions could lead to exposure to infectious agents.

**DANGER**

Reusable accessories should be cleaned and disinfected or sterilized as stated by the manufacturer, after each patient examination.

**WARNING**

- Be sure to verify the media after writing data, such as EZBackup, SaveAs or Export.
- Before deleting a patient or image from the patient screen, make sure you have saved the data by EZBackup/Backup or Export and verify that the media data transfer was successful. Not following the instructions could result in the loss of data requiring a rescan.

**WARNING**

When you move the Control Panel up/down with the monitor, place BOTH hands on the Control Panel. Touching other moving parts other than the Control Panel may cause personal injury.

**WARNING**

Contact with natural rubber latex may cause a severe anaphylactic reaction in persons sensitive to the natural latex protein. Sensitive users and patients must avoid contact with these items.

**WARNING**

The ultrasound system is not intended to be used as a data storage device; backup of the Patient and Image Database is your institution's responsibility. GE HealthCare is NOT responsible for any lost patient information or for lost images. Loss of image data may require a rescan.

**WARNING**

To minimize accidental loss of data, perform EZBackup and Backup on a regular basis.

1. First, perform EZBackup to save the images.
2. Next, perform Backup at **Utility > Backup/Restore**. Enable the following checkboxes under Backup:
 - User defined configuration
 - Service

**CAUTION**

Only approved and recommended peripherals and accessories should be used.

All peripherals and accessories must be securely mounted to the ultrasound system. Failure to follow these instructions could lead to unexpected diagnostic performance.

**CAUTION**

Non-supported peripheral devices that use their own AC power source CANNOT be attached to the ultrasound system. DO NOT connect the peripheral device's power cord into the ultrasound system. Only peripheral devices purchased from GE HealthCare with the purpose of being used with the ultrasound system should be used.

Use a USB printer cable that is less than 3 meters in length.

Failure to follow these instructions could lead to unexpected diagnostic performance.

**CAUTION**

Do not use this equipment if a safety problem is known to exist. Have the unit repaired and performance verified by qualified service personnel before returning to use.

**CAUTION**

To avoid injury or system damage, NEVER place any object or liquid on the operator panel.

**CAUTION**

Archived data is managed at the individual sites. Performing data backup (to any device) is recommended.

**CAUTION**

DO NOT touch the patient and any of the connectors on the ultrasound unit simultaneously, including ultrasound probe connectors.

DO NOT touch the conducting parts of the USB, Ethernet, Video, HDMI, Audio cables when connecting equipment to the unit.

Material Safe Data

Rubber Part

Material: GFE 15&17

Where Used: Probe holder/Rubber cover for arm stopper/ Rubber plug

Monitor Related Hazards



WARNING

- DO NOT place a finger, hand or any object on the joint of the monitor or monitor arm to avoid injury when moving the monitor and monitor arm.
- To avoid result of injury or system damage, NEVER place any object or liquid on the monitor, whether in the home or flip down/transport position.
- DO NOT scratch or press on the panel with any sharp objects, such as a pencil or pen, as this may result in damage to the panel.
- To avoid injury or damage, make sure nothing is within the range of motion before moving the monitor and monitor arm. This includes both objects and people.
- Pay attention to the monitor arm position to avoid hitting it against anyone or anything.
- Before moving the system to another location, be sure to lock the monitor arm in the transport position.
- The monitor screen may have defective pixels. These pixels may appear as a slightly light or dark area on the screen. This is due to the characteristics of the panel itself, and not the product.
- The backlight of the monitor has a fixed life span. When the screen becomes dark or begins to flicker, contact a qualified Service Representative for information.

NOTE

Bright light could impact readability of screen.

Material Safe Data

Rubber Part

Material: Silicon

Where Used: Handle Screw Cap/Monitor Rubber

Allergic reactions to latex-containing medical devices



CAUTION

Due to reports of severe allergic reactions to medical devices containing latex (natural rubber), the FDA advises health-care professionals to identify latex-sensitive patients, and be prepared to treat allergic reactions promptly. Latex is a component of many medical devices, including surgical and examination gloves, catheters, incubation tubes, anesthesia masks and dental dams. Patient reaction to latex has ranged from contact urticaria, to systemic anaphylaxis.

For more details regarding allergic reaction to latex, refer to FDA Medical Alert MDA91-1, March 29.

Electrical safety

Device classifications

Equipment is Class I or Internally Powered equipment, Type BF or Type CF Applied Part for probes, DEFIBRILLATION-PROOF Type CF Applied Part for ECG.

Internally connected peripheral devices

The system, together with peripheral devices, such as video printer, meets IEC 60601-1 standards for electrical isolation and safety. These standards are applicable only when the specified peripheral devices are plugged into the AC outlets provided on the unit.

External Connection of other peripheral devices



CAUTION

External devices can be used only if CE marked and in compliance with related standards (EN 60601-1 or EN 60950/62368). Conformance to EN 60601-1 must be verified.

Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards (e.g. IEC 60950/62368 for data processing equipment and IEC 60601-1 for medical equipment). Furthermore all complete configurations shall comply with the valid version of the system standard IEC 60601-1. Anybody connecting additional equipment to the signal input part or signal output part of the ultrasound unit configures a medical system, and is therefore responsible that the system complies with the requirements of the valid version of IEC 60601-1. If in doubt consult the technical service department or your local GE HealthCare representative.

Other external devices, such as printers, and external monitors, usually exceed allowable leakage limits and, when plugged into separate AC outlets that are then connected to the unit, are in violation of patient safety standards. Suitable electrical isolation of such external AC outlets may be required in order to meet UL60601-1 and IEC 60601-1 standards for electrical leakage.

EMC (Electromagnetic Compatibility)

The product was tested according to the recommendations of IEC TR 60601-4-2: Medical electrical equipment – Part 4-2: Guidance and interpretation – Electromagnetic immunity: performance of medical electrical equipment and medical electrical systems.

NOTE

This equipment generates, uses and can radiate radio frequency energy. The equipment may cause radio frequency interference to other medical and non-medical devices and radio communications. To provide reasonable protection against such interference, this product complies with emissions limits for a Group 1, Class A Medical Devices Directive as stated in EN 60601-1-2. However, there is no guarantee that interference will not occur in a particular installation.

NOTE

If this equipment is found to cause interference (which may be determined by turning the equipment on and off), the user (or qualified service personnel) should attempt to correct the problem by one or more of the following measure(s):

- Reorient or relocate the affected device(s)
- Increase the separation between the equipment and the affected device
- Power the equipment from a source different from that of the affected device
- Consult the point of purchase or service representative for further suggestions.

NOTE

The manufacturer is not responsible for any interference caused by using other than recommended interconnect cables or by unauthorized changes or modifications to this equipment. Unauthorized changes or modifications could void the users' authority to operate the equipment.

NOTE

To comply with the regulations on electromagnetic interference for a Class A FCC Device, all interconnect cables to peripheral devices must be shielded and properly grounded. Use of cables not properly shielded and grounded may result in the equipment causing radio frequency interference in violation of the FCC regulations.

"Harmful interference" is defined in 47 CFR §2.122 by the FCC as follows: Interference which endangers the functioning of a radionavigation service or of other safety services or seriously degrades, obstructs, or repeatedly interrupts a radio communication service operating in accordance with the [ITU] Radio Regulations.

NOTE

Do not use devices which intentionally transmit RF Signals (cellular phones, transceivers, or radio controlled products), other than those supplied by GE HealthCare, in the vicinity of the equipment, as it may cause performance outside the published specifications. Keep the power to these type devices turned off when near this equipment.

The medical staff in charge of this equipment is required to instruct technicians, patients, and other people who maybe around this equipment to fully comply with the above requirement.

EMC Performance

All types of electronic equipment may characteristically cause electromagnetic interference with other equipment, either transmitted through air or connecting cables. The term EMC (Electromagnetic Compatibility) indicates the capability of equipment to curb electromagnetic influence from other equipment and at the same time not affect other equipment with similar electromagnetic radiation from itself.

Proper installation following the service manual is required in order to achieve the full EMC performance of the product.

The product must be installed as stipulated in Notice upon Installation of Product section.



CAUTION

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally. Failure to follow these instruction could lead to unexpected diagnostic performance.

In case of issues related to EMC, please call your service personnel.

The manufacturer is not responsible for any interference caused by using other than recommended interconnect cables or by unauthorized changes or modifications to this equipment. Unauthorized changes or modifications could void the users' authority to operate the equipment.



WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ultrasound system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.



WARNING

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.



WARNING

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Cable Specification

This device contains the following cables:

Table 4-2 Cable specification list

Number	Name	Cable Length	Cable Shielded
1	Power supply cord	1m<L<3m	Unshielded
2	Probe cord (except for 6Tc-RS)	1m<L<3m	Shielded
3	6Tc-RS probe cord	>3m	Shielded

Number	Name	Cable Length	Cable Shielded
4	HDMI cable for connection to Digital Expert Connect, external display etc.	1m<L<3m	Shielded
5	HDMI cable option for connection to Video Output Adapter	0.375m<L <1m	Shielded
6	Type-C cable for connection to external device	1m<L<3m	Shielded
7	Ethernet cable for connection to local area network	1m<L<3m	Unshielded
8	Port on Video adapter, connect through Video adapter HDMI transmission cable	1m<L<3m	Shielded
9	Port on Video adapter, connect through S-video transmission cable	3m	Shielded
10	Port on Video adapter, connect through CVBS transmission cable	3m	Shielded
11	Power supply cord for Gel Warmer	<1m	Unshielded
12	NORAV ECG lead wire	>3m	Shielded
13	NORAV ECG USB Cable	0.375m<L <1m	Shielded
14	Footswitch MKF 2 1S/1S-MED USB GP26 transmission cable	1m<L<3m	Shielded
15	Footswitch FSU-1000 transmission cable	1m<L<3m	Shielded
16	Western Digital 1TB mobile USB HDD transmission cable	0.375m<L <1m	Shielded
17	Transcend TS8XDVDS-K transmission cable	0.375m<L <1m	Shielded
18	NETGEAR wireless adapter A8000 transmission cable	0.375m<L <1m	Shielded
19	SONY UP-D25MD Printer transmission cable	1m<L<3m	Shielded
20	SONY UP-D898MD Printer transmission cable	1m<L<3m	Shielded
21	Printer OfficeJet 200 Mobile Printer transmission cable	1m<L<3m	Shielded
22	Printer USB isolator transmission cable	0.375m<L <1m	Shielded
23	Barcode Reader 1950 transmission cable	1m<L<3m	Shielded

Notice upon Installation of Product

Separation distance and effect from fixed radio communications equipment: field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and

land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast transmitter cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the ultrasound system is used exceeds the applicable RF compliance level as stated in the immunity declaration, the ultrasound system should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as re-orienting or relocating the ultrasound system or using an RF shielded examination room may be necessary.

1. Use either power supply cords provided by GE HealthCare or ones designated by GE HealthCare. Products equipped with a power source plug should be plugged into the fixed power socket which has the protective grounding conductor. Never use any adaptor or converter to connect with a power source plug (e.g. three-prong-to-two-prong converter).
2. Locate the equipment as far away as possible from other electronic equipment.
3. Be sure to use only the cables provided by or designated by GE HealthCare. Connect these cables following the installation procedures (e.g. wire power cables separately from signal cables).
4. Lay out the main equipment and other peripherals following the installation procedures described in the Option Installation manuals.

General Notice



WARNING

Use of accessories, transducers and cables other than those specified or provided by GE of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

1. Designation of Peripheral Equipment Connectable to This Product.
 The equipment indicated in the Supplies/Accessories section can be hooked up to the product without compromising its EMC performance.
 Avoid using equipment not designated in the list. Failure to comply with this instruction may result in poor EMC performance of the product.
2. Notice against User Modification.
 The user should never modify this product. User modifications may cause degradation in EMC performance.
 Modification of the product includes changes in:
 Cables (length, material, wiring, etc.)
 System installation/layout
 System configuration/components
 Securing system parts (cover open/close, cover screwing)
3. Operate the system with all covers closed. If a cover is opened for some reason, be sure to shut it before starting/resuming operation.
4. Operating the system with any cover open may affect EMC performance.

Peripheral Update for EC countries

The following is intended to provide the users in EC countries with updated information concerning the connection of the ultrasound system to image recording and other devices or communication networks.

Peripherals used in the patient environment

The ultrasound system has been verified for overall safety, compatibility and compliance with the image recording devices listed in Supplies/Accessories section.

The ultrasound system has also been verified for compatibility, and compliance for connection to a local area network (LAN) via the rear panel Ethernet connection, provided the LAN components are IEC/EN 60950/62368 compliant.

The ultrasound system may also be used safely while connected to devices other than those recommended above if the devices and their specifications, installation, and interconnection with the system conform to the requirements of IEC/EN 60601-1.

Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards (i.e., IEC 60950/62368 for data processing equipment and IEC 60601-1 for medical equipment). Furthermore, all complete configurations shall comply with the valid version of the system standard IEC 60601-1. Everyone who connects additional equipment to the signal input part or signal output part of the ultrasound system configures a medical system, and is therefore responsible to ensure that the system complies with the requirement of the valid version of IEC 60601-1. If in doubt, consult the technical service department or your local GE representative.

General precautions for installing an alternate on-board device would include:

1. The added device(s) must have appropriate safety standard conformance and CE Marking.
2. The total power consumption of the added devices, which connect to the ultrasound system and are used simultaneously, must be less than or equal to the rated supply of the ultrasound system.
3. There must be adequate heat dissipation and ventilation to prevent overheating of the device.
4. There must be adequate mechanical mounting of the device and stability of the combination.
5. Risk and leakage current of the combination must comply with IEC/EN 60601-1.
6. Electromagnetic emissions and immunity of the combination must conform to IEC/EN 60601-1-2.
7. The added device(s) must be used for their intended purpose having a compatible interface.
8. Signal or mains isolation devices and additional protective earth may be needed to assure compliance with IEC/EN 60601-1.



CAUTION

The connection of equipment or transmission networks other than as specified in the user instructions can result in an electric shock hazard or equipment malfunction. Substitute or alternate equipment and connections requires verification of compatibility and conformity to IEC/EN 60601-1 by the installer. Equipment modifications and possible resulting malfunctions and electromagnetic interference are the responsibility of the owner.

Peripheral used in the non-patient environment

The ultrasound system has also been verified for compatibility, and compliance for connection to a USB HDD/USB memory via the system USB port, provided the USB HDD/USB memory are IEC/ EN 60950/62368 compliant.

Interference caution

Intended healthcare environment: Professional healthcare facility environment.



CAUTION

Use of devices that transmit radio waves near the unit could cause it to malfunction.

Devices which intrinsically transmit radio waves, such as cellular phones, radio transceivers, mobile radio transmitters, radio-controlled toys, etc., should preferably not be operated near the unit. See *Minimum distances* on page 278 about the recommended minimum separation distances between portable and mobile RF communications equipment and the ultrasound unit.

Medical staff in charge of the unit are required to instruct technicians, patients, and other people who may be around the unit to fully comply with the above recommendations.

Any electrical device can unintentionally emit electromagnetic waves. However, minimum device separation distances cannot be calculated for such unspecified radiation. When the ultrasound unit is used adjacent to or in close proximity to other equipment the user should be attentive to unexpected device behavior which may be caused by such radiation.

The Medical Electrical System is suitable for use in professional healthcare facility environment.

The ultrasound unit is intended for use in the electromagnetic environment specified in the tables below.

The user of the ultrasound unit should ensure that the device is used in such an environment.

Declaration of Emissions

This system is suitable for use in the following environment. The user must assure that it is used only in the electromagnetic environment as specified.

Table 4-3 Declaration of Emissions

Guidance and manufacturer's declaration - electromagnetic emissions		
The system is intended for use in the electromagnetic environment specified below. The user of the system should assure that it is used in such an environment.		
Emission Type	Compliance	Electromagnetic Environment
RF Emissions CISPR 11	Group 1	This system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

Guidance and manufacturer's declaration - electromagnetic emissions		
The system is intended for use in the electromagnetic environment specified below. The user of the system should assure that it is used in such an environment.		
Emission Type	Compliance	Electromagnetic Environment
RF Emissions CISPR 11	Class A	This system is suitable for use in all establishments, other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	

**WARNING**

The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the equipment.

**WARNING**

This system is intended for use by healthcare professionals only. This system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the system or shielding the location.

Declaration of Immunity

The system is suitable for use in the specified electromagnetic environment and it has meets the following immunity test levels. Higher immunity levels may cause the system's essential performance lost or degraded.

Table 4-4 Electromagnetic Immunity

Phenomenon	Basic EMC Standard or Test Method	Regulatory Acceptable Level	EMC Environment and Guidance
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. Mains power quality should be that of a typical commercial and/or hospital environment. If the user requires continued operation during power mains interruptions, it is recommended that the system be powered from a UPS or a battery. NOTE: UT is the a.c. mains voltage prior to application of the test level. Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial and/or hospital environment.
Radiated RF EM fields	IEC 61000-4-3	3 V/m 80 MHz - 2.7GHz 80%AM at 1kHz	
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See the RF wireless communication equipment table in "Minimum distances".	
Rated power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	
Electric fast transients bursts	IEC 61000-4-4	Power supply port: ± 2 kV, 100 kHz repetition frequency; Signal input/output parts Port: ± 1 kV, 100 kHz repetition frequency.	

Table 4-5 Electromagnetic Immunity (continued)

Phenomenon	Basic EMC Standard or Test Method	Regulatory Acceptable Level	EMC Environment and Guidance
Surges	IEC 61000-4-5	Line to line: ± 0.5 kV, ± 1 kV Line to earth: ± 0.5 kV, ± 1 kV, ± 2 kV	Separation distance to radio communication equipment must be maintained according to the method below. Interference may occur in the vicinity of equipment marked with the symbol:  Image degradation or interference may occur due to conducted RF noise on the equipment mains power supply or other signal cable. Such interference is easily recognized and distinguishable from patient anatomy and physiological waveforms. Interference of this type may delay the examination without affecting diagnostic accuracy. Additional mains/signal RF isolation or filtering may be needed if this type interference occurs frequently.
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V in 0.15 MHz - 80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz	
Voltage dips	IEC 61000-4-11	0% UT: 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% UT: 1 cycle and 70% UT: 25/30 cycles sine phase at 0°	
Voltage interruptions	IEC 61000-4-11	0% UT: 250/300 cycle	
NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people. If noise generated from other electronic equipment is near the probe's center frequency, noise may appear on the image. Good power line isolation is required.			

Minimum distances

Recommended separation distances between portable and mobile RF communications equipment and the ultrasound unit

The ultrasound unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the ultrasound unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ultrasound unit as recommended below, according to the maximum output power of the communications equipment. Portable and mobile radio communications equipment (e.g. two-way radio, cellular/cordless telephones and similar equipment) should be used no closer to any part of this system, including cables, than determined according to the following method:

Recommended separation distances between portable and mobile RF communications equipment and the ultrasound unit						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380 - 390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	430 - 470	GMRS 460 FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 -787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 - 960	GSM 800/900, tetra 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700 - 1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25, UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100 - 5800	WLAN, 802.11 a/n,	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						
For transmitters rated at a maximum output power not listed above the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.						

Wireless Adapter Specification

Table 4-6 Wireless Adapter specification

Received frequency (GHz)	Emission frequency (GHz)	Modulation Type	Effective radiation frequency
2.412	2.412	OFDM	<=100mW
2.417	2.417		
2.422	2.422		
2.427	2.427		
2.432	2.432		
2.437	2.437		
2.442	2.442		
2.447	2.447		
2.452	2.452		
2.457	2.457		
2.462	2.462		

Bluetooth Specification

- Frequency Range: 2.4GHz ISM Band (2.402 - 2.480 GHz)
- Maximum Output Power: ≤100 mW (+20 dBm)
- Data Transmission Rate: 3 M/s (10 meters on the PC to PC)
- Modulation Type: FSK, GFSK
- Transmission protocol: Bluetooth protocol



WARNING

Even if other devices meet the emission requirements of the corresponding national standards, the device or system may still be interfered by other devices.

Essential performance

The essential performance of the ultrasound unit is:

1. The ability to be free from noise on a waveform or artefacts or distortion in an image or error of a displayed numerical value which cannot be attributed to a physiological effect and which may alter the diagnosis.
2. The ability to be free from the display of incorrect numerical values associated with the diagnosis to be performed ^a.
3. The ability to be free from the display of incorrect safety-related indications ^a.
4. The ability to be free from the production of unintended or excessive ultrasound output.
5. The ability to be free from the production of unintended or excessive TRANSDUCER ASSEMBLY surface temperature.
6. The ability to be free from the production of unintended or uncontrolled motion of TRANSDUCER ASSEMBLIES intended for intra-corporeal use.

NOTE

^a: “incorrect” in the sense that the displayed value differs from what is calculated (having been altered during data transfer), or the calculation itself is not correct.

Performance as per IEC TR 60601-4-2

Performance as per IEC TR 60601-4-2 for the ultrasound unit is:

- The ability to display physiological images as input for diagnosis by trained physician.
- The ability to display physiological traces as aid for diagnosis by trained physician.
- The ability to display quantified data as input for diagnosis by trained physician.
- No reset to factory defaults (manufacturer’s presets).

Patient Environmental Devices

1. Left side
 - Gel bottle holder / Gel warmer
 - 2 probe holders
2. Right side
 - Probe holders
3. Front panel
 - Inside Printer Box

- 4 or 5 active probe ports (For Versana Premier Lotus)
 - 3 or 4 active probe ports (For Versana Premier)
4. Rear panel
- 1 HDMI port (1080P)
 - 3 USB ports (2 type A ports, 1 type C port)
 - 1 Ethernet port (100MBase-T)

Acceptable Devices

The Patient Environmental devices shown on the previous page are specified to be suitable for use within the PATIENT ENVIRONMENT.



CAUTION

DO NOT connect any probes or accessories without approval by GE HealthCare within the PATIENT ENVIRONMENT.

Unapproved Devices



CAUTION

DO NOT use unapproved devices.

If devices are connected without the approval of GE HealthCare, the warranty will be INVALID.

Any device connected to the Versana Premier/Versana Premier Lotus must conform to the requirements for IEC or equivalent standards appropriate to devices.

Accessories, Options, Supplies



CAUTION

Unsafe operation or malfunction may result. Use only the accessories, options and supplies approved or recommended in these instructions for use.

Acoustic Output



CAUTION

Allowing the machine to transmit acoustic output with the probe not in use (or in its holder) can cause the transducer to build up heat. Always lower the acoustic power or freeze the image when not in use.

When the **Auto Freeze** preset is selected on the **Utility > System > System Imaging** screen, the system auto freezes if it detects no change in the image.

Located on the upper right section of the system display monitor, the acoustic output display provides the operator with real-time indication of acoustic levels being generated by the system, and the acoustic output level could be adjusted by user. See the Acoustic Output chapter in the Advanced Reference Manual for more information.

Acoustic Output Display Specifications

The display consists of three parts: Thermal Index (TI), Mechanical Index (MI), and a relative Acoustic Output (AO) value. Although not part of the IEC standard, the AO value informs the user of where the system is operating within the range of available output.

The TI and MI are displayed at all times. The TI and MI display starts at a value of 0 and increments in steps of ≤ 0.2 .

Always be aware of the acoustic output level by observing the Acoustic Output Display. In addition, become thoroughly familiar with the Acoustic Output Display and equipment controls affecting output.

Thermal Index

Depending on the examination and type of tissue involved, the TI parameter will be one of three types:

- Soft Tissue Thermal Index (TIS). Used when imaging soft tissue only, it provides an estimate of potential temperature increase in soft tissue.
- Bone Thermal Index (TIB). Used when bone is near the focus of the image as in the third trimester OB examination, it provides an estimate of potential temperature increase in the bone or adjacent soft tissue.
- Cranial Bone Thermal Index (TIC). Used when bone is near the skin surface as in transcranial examination, it provides an estimate of potential temperature increase in the bone or adjacent soft tissue.

Mechanical Index

MI recognizes the importance of non-thermal processes, cavitation in particular. The Index is a relative indicator of the likelihood of mechanical bioeffect within the tissue.

Changing the Thermal Index Type

You can select the displayed TI type on Utility > Imaging > B-Mode. This preset is application dependent so each application could specify a different TI type.

TI and MI Display Accuracy

When display $MI \geq 0.6$, $TI \geq 3.6$, the displayed values of MI and IT is not lower than 50% or higher than 150% of the measured value.

When display $MI < 0.6$, $TI < 3.6$, the absolute error of $MI \leq 0.3$, the absolute error of $TI \leq 1.8$.

Controls Affecting Acoustic Output

The potential for producing mechanical bioeffects (MI) or thermal bioeffects (TI) can be influenced by certain controls.

Direct. The Power Output control has the most significant effect on Acoustic Output.

Indirect. Indirect effects may occur when adjusting controls. Controls that can influence MI and TI are detailed under the Bioeffects portion of each control in the Optimizing the Image sections.

Always observe the Acoustic Output display for possible effects.

Best Practices While Scanning



HINTS

Raise the Acoustic Output only after attempting image optimization with controls that have no effect on Acoustic Output, such as Gain and TGC.

NOTE

Refer to the Optimizing the Image sections for a complete discussion of each control.



DANGER

Be sure to have read and understood control explanations for each mode used before attempting to adjust the Acoustic Output control or any control that can affect Acoustic Output. During a screening and diagnostic ultrasound examination, high frequency sound penetrates and interacts with tissue in and around the area of anatomy to be imaged. Only a small portion of this sound energy is reflected back to the transducer for use in constructing the image, while the remainder is dissipated within the tissue. The interaction of sound energy with tissue at sufficiently high levels can produce biological effects (aka bioeffects) of either a mechanical or thermal nature. Although the generation of bioeffect is intentional with therapeutic ultrasound, it is generally undesired in screening and diagnostic applications and may be harmful in some conditions.



ACOUSTIC OUTPUT HAZARD

Use the minimum necessary acoustic output to get the best diagnostic image or measurement during an examination. Begin the exam with the probe that provides an optimum focal depth and penetration.

Acoustic Output Default Levels

In order to assure that an exam does not start at a high output level, the ultrasound system initiates scanning at a reduced default output level. This reduced level is preset programmable and depends upon the exam category and probe selected. It takes effect when the system is powered on or Patient is selected.

To modify acoustic output, adjust the Power Output level.

Safety statement

Although no harmful biological effects have been demonstrated for ultrasound frequencies, intensities, or exposure times used in examination with the GE HealthCare system, GE HealthCare recommends using the lowest acoustic output settings which will produce diagnostically acceptable information.

RoHS Hazardous Substances

The following product pollution control information is provided according to SJ/T11364-2014 Marking for Control of Pollution caused by Electronic Information Products.



The symbol above indicates the product contains hazardous materials in excess of the limits established by the Chinese standard GB/T 26572 Requirements for Concentration Limits for Certain Hazardous Substances in Electronic Information Products. The number in the symbol is the Environment-friendly Use Period (EFUP), which indicates the period during which the hazardous substances or elements contained in electrical and electronic products will not leak or mutate under normal operating conditions so that the use of such electrical and electronic products will not result in any severe environmental pollution, any bodily injury or damage to any assets. The unit of the period is "Year."

In order to maintain the declared EFUP, the product shall be operated normally according to the instructions and environmental conditions as defined in the product manual, and periodic maintenance schedules specified in Product Maintenance Procedures shall be followed strictly. Consumables or certain parts may have their own label with an EFUP value less than the product. Periodic replacement of those consumables or parts to maintain the declared EFUP shall be done in accordance with the Product Maintenance Procedures.

This product must not be disposed of as unsorted municipal waste, and must be collected separately and handled properly after decommissioning.

Name and Concentration of Hazardous Substances

Table 4-7 Table of hazardous substances' name and concentration

Component Name	Hazardous substances' name					
	Pb	Hg	Cd	Cr (VI)	PBB	PBDE
LCD Monitor	X	O	O	O	O	O
Printed Circuit Board Assemblies	X	O	O	O	O	O
Touch Panel	X	O	O	O	O	O
Keyboard Assemblies	X	O	O	O	O	O
Power Assemblies	X	O	O	O	O	O
Console Cabinet	O	O	O	O	O	O
Ultrasound Probes	X	O	O	O	O	O
Wheels	O	O	O	O	O	O

This table is prepared according to SJ/T 11364.

O: Indicates that this hazardous substance contained in all of the homogeneous materials for this part is below the limit requirement in GB/T 26572.

X: Indicates that this hazardous substance contained in at least one of the homogeneous materials used for this part is above the limit requirement in GB/T 26572.

Data listed in the table represents best information available at the time of publication.

Applications of hazardous substances in this medical device are required to achieve its intended clinical uses, and/or to provide better protection to human beings and/or to environment, due to lack of reasonably (economically or technically) available substitutes

Note: Options may not be present on every system.

WEEE Passport

The WEEE (Waste Electrical and Electronic Equipment) Passport describes product recycling information. To access the WEEE passport for GE HealthCare products:

1. Go to the GE HealthCare Support Documentation Library at:
<https://www.gehealthcare.com/support/documentation>
2. Select the modality "Ultrasound (UL)."
3. Enter the document name or the keyword "WEEE."
4. Press "Search."
5. Select the desired WEEE passport.

Safe Product and Packaging Disposal

This product and package should be disposed of according to hospital disposal practices, and local environmental and waste disposal regulations. Components and accessories of the ultrasound system which have come into direct or indirect contact with the patient may be biohazardous, and should be disposed of according to facility guidelines for biohazardous material. The waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the disposal/decommissioning of equipment.

General Caution



CAUTION

Standard maintenance must be performed by authorized service personnel for the lifetime of the product (7 years).



CAUTION

Proceed cautiously when crossing door or elevator thresholds. Use the handle to push/pull the system, e.g., do not use the monitor. Failure to do so may cause serious injury or system damage.

Device Labels

Label Icon Description

The following table describes the purpose and location of safety labels and other important information provided on the equipment.

NOTE

This machine should be used in compliance with law. Some jurisdictions restrict certain uses, such as gender determination (shown for country specific label).

Table 4-8 Label Icons

Label/Icon	Purpose/Meaning	Location	Reference
	Manufacture's name and address	Rating Plate	EN ISO 15223-1
	Date of manufacture: The date could be a year, year and month, or year, month and day, as appropriate.	Rating Plate	EN ISO 15223-1
	Serial Number	Rating Plate	EN ISO 15223-1
	Catalog Number	Rating Plate	EN ISO 15223-1
	Part Number	Rating Plate	GE HealthCare created
IP Code: IPX8 FSU-1000, MKF 2-MED GP26	Indicates the degree of protection provided by the enclosure per IEC 60529. Protected against the effects of continuous immersion.	Bottom of footswitch	IEC 60601-1 (EN 60601-1)
IP Code: IPX7 Probe head	Indicates the degree of protection provided by the enclosure per IEC 60529. Protected against the effects of temporary immersion.	Probe Connector	IEC 60601-1 (EN 60601-1)
	Authorized European Representative address.	Rating Plate	EN ISO 15223-1

Label/Icon	Purpose/Meaning	Location	Reference
	The CH-REP label indicates the authorized representative in Switzerland.	Rating Plate	Swiss Medical Devices Ordinance (MedDO)
	Type CF Defib-Proof Applied Part (heart in the box with paddle) symbol is in accordance with IEC 60417-5336.	ECG Bracket and ECG Module	IEC 60601-1 (EN 60601-1)
	General Warning.	Various	IEC 60601-1 (EN 60601-1)
	Avoid pinching hands during adjusting the monitor or control panel height adjustment assy.	Various	GE HealthCare created
	"CAUTION" - Dangerous voltage" (the lightning flash with arrowhead) is used to indicate electric shock hazards.	Body back cover	IEC 60601-1 (EN 60601-1)
	Indicates the power on (green), power off (white) and no power (no light) position. CAUTION: This Power Switch DOES NOT ISOLATE Mains Supply.	See the Console Overview section for location information	IEC 60601-1 (EN 60601-1)
	"Protective Earth" indicates the protective earth (grounding) terminal.	Inside Power Box	IEC 60601-1 (EN 60601-1)
	"Equipotentiality" indicates the terminal to be used for connecting equipotential conductors when interconnecting (grounding) with other equipment.	Rear panel	IEC 60601-1 (EN 60601-1)
	This symbol indicates that waste electrical and electronic equipment must not be disposed of as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.	Body back cover and probe connector	WEEE Recast Directive 2012/19/EU

Label/Icon	Purpose/Meaning	Location	Reference
	The separate collection symbol is affixed to a battery, or its packaging, to advise you that the battery must be recycled or disposed of in accordance with local or country laws. Please consult the service manual or equipment instructions. Information on the potential effects on the environment and human health of the substances used in batteries is available at this url: http://www.gehealthcare.com/euen/wEEE-recycling/index.html	Battery Pack	2006/66/EC Directive
	Do not place any objects on the monitor.	Back of Monitor	GE HealthCare created
	Do not force the monitor with your hands.	Back of Monitor	GE HealthCare created
	There is a pinch point on the monitor. Take care to avoid injuring hands or fingers when flipping down the monitor.	Back of Monitor	GE HealthCare created
	Use the rear handle for horizontal movement only.	Back of Monitor	GE HealthCare created
	"Consult accompanying documents" is intended to alert the user to refer to the operator manual or other instructions when complete information cannot be provided on the label.	Various	EN 60601-1/A1 (IEC 60601-1/A1)

Label/Icon	Purpose/Meaning	Location	Reference
	DO NOT place a finger, hand or any object on the joint of the monitor or monitor arm to avoid injury when moving the monitor and monitor arm.	Flexible Arm	GE HealthCare created
	The maximum load weight is 1.5kg.	Various	GE HealthCare created
	Do not push the system body without using handle.	Back of Monitor	EN 60601-1/A1 (IEC 60601-1/A1)
	GOST symbol: Russia Regulatory Country Clearance. Note: Only after Russian regulatory registration is complete, this label will be located on the console rating plate.	Rating plate	Russia regulatory requirements.
	Ukrainian mark of conformity with the Technical Regulations. This product meets the requirements of the Technical Regulations on medical devices, approved by Resolution No.753 of the Cabinet of Ministers of Ukraine of 2 October 2013. Note: Only after Ukrainian regulatory registration is complete, this label will be located on the console.	Body back cover	Ukraine regulatory requirements.
	Indian BIS certificate Mark. Note: Only after Indian BIS certificate is complete, this icon will be located on the battery or LCD.	Battery	Indian BIS requirements.
	INMETRO Certification: SGS Note: Only after Brazilian regulatory registration is complete, this label will be located on the console rating plate.	Rating plate	Brazil regulatory requirement.
	CAUTION: This machine should be used in compliance with law. Some Jurisdictions restrict certain uses, such as gender determination. Note: For China, Korea and India only.	Body back cover	GE HealthCare created

Label/Icon	Purpose/Meaning	Location	Reference
	The maximum weight of the system with all hardware and mechanical configurations connected.	Rating plate	GE HealthCare created
	Every system has a unique marking for identification, the Unique Device Identification (UDI) Label. The UDI label consists of a series of alpha-numeric characters and barcode which uniquely identify the system as a medical device manufactured by General Electric. Scan or enter the UDI information into the patient health record as required by country-specific laws.	Rating plate	GE HealthCare created
	Federal law restricts this device to sale by or on the order of a physician.	Body back cover	United States regulatory requirement.
	The CE Mark of Conformity indicates this equipment conforms with the MDR - REGULATION (EU) 2017/745.	Rating plate	Europe regulatory requirement.
	Input	Rating plate	GE HealthCare created
	This symbol indicates that the device is a Medical Device according to MDR - REGULATION (EU) 2017/745.	Rating plate	EN ISO 15223-1
	The symbol indicates that the device is only for GE HealthCare Ultrasound device.	Gel Warmer	GE HealthCare created
	Symbol indicating that the Instructions for Use are supplied in electronic form. NOTE: The symbol is valid for eIFU only when the system meets eIFU requirement.	Body back cover	ISO7000-3500
	Keep the system away from magnetic resonance imaging (MRI) equipment.	Body back cover	ASTM F2503-20

Label Locations

Ultrasound system warning labels are provided in English.




Table 4-9 Label Location Explanations

No	Explanation
1	Possible shock hazard. Do not remove covers or panels. No user serviceable parts are inside. Refer servicing to qualified service personnel.
2	MRI Warning label.
3	Do not use the following devices near this equipment: cellular phone, radio receiver, mobile radio transmitter, radio controlled toy, etc. Use of these devices near this equipment could cause this equipment to perform outside the published specifications. Keep power to these devices turned off when near this equipment.
4	WEEE Label
5	United States only Prescription Requirement label
6	Symbol indicating that the Instructions for Use are supplied in electronic form.
7	"Consult accompany document" is intended to alert the user to refer to the operator manual or other instructions when complete information cannot be provided on the label.
8	NRTL Listing and Certification Mark.

No	Explanation
9	CISPR CAUTION: The ultrasound system conforms to the CISPR11, Group 1, Class A of the international standard for Electromagnetic disturbance characteristics.

NOTE

The warning labels are on the rear of the system. The label content may be different for different regions. Please refer to the warning labels on the system for the actual content.

Figure 4-1 Label Location - example



Figure 4-2 Battery Label



1. Do not put the battery in fire.
2. Do not disassemble or mistreat the battery.



CAUTION

Do not disassemble or mistreat the battery. Do not put the battery in fire. Replace the battery with the same battery type only. Failure to follow these instructions may present risk of explosion fire or high temperature. See the battery user manual for additional safety instructions.

Probe Label Explanation

The following information appears on all probe labels, regardless of the connector type, except for "IPX7", "CE Mark", and "XDclear™" which only appears on applicable probes.

NOTE

The probe label displayed in this manual is only for illustrational purposes. Refer to the actual probe label for the specific information.

Figure 4-3 Probe label




Table 4-10 Probe label icon explanation

1. GE Logo 2. Probe Model (Name) 3. UDI Symbol and Data Matrix 4. UDI Human Readable Label Text: Global Trade Item Number, GTIN, (01), Manufacturing Date (11), Serial Number (21) 5. Probe Applied Part (Can be CF or BF Type: depends on applied part symbol used.) 6. Caution: Consult the Manual 7. WEEE Waste Symbol	8. Chinese RoHS Hazardous Substance Symbol 9. CE Mark and Notified Body Number 10. REF: Catalog/model number 11. Serial Number 12. Manufacturer's site country of origin 13. Legal Manufacturer's Name and Address 14. Date of Manufacture, as YYYY-MM 15. IP Classification
--	--

NOTE

Non-GE probes will also have a UDI symbol and equivalent information.

Probe Box Label

NOTE

The probe box label displayed in this manual is only for illustrational purposes. Refer to the actual probe box label for the specific information.

Figure 4-4 Probe box label

Cat#		 * H 4 5 0 2 1 R P *		YEAR 16 QTR 2 CE 0459
H45021RP				
Desc: 4C-RS		Service Part#: 5499316  * 5 4 9 9 3 1 6 *		
SN  * 1 2 3 4 5 6 P D 1 *		123456PD1		
UDI  (01)00840682111157(11)160400(21)123456PD1				

1. Probe Box Label UDI GTIN
2. Probe Box Label Barcode Location

Chapter 5

Probes and Biopsy

In this section

Probe Overview.	298
Probe Discussion.	329
Biopsy Special Concerns.	339
Preparing for a Biopsy.	341

Probe Overview

Ergonomics

Probes have been ergonomically designed to:

- Handle and manipulate with ease
- Connect to the system with one hand
- Be lightweight and balanced
- Have rounded edges and smooth surfaces
- Stand up to typical wear by cleaning and disinfectant agents, contact with approved gel, etc.

Cables have been designed to:

- Connect to system with appropriate cable length

Cable handling

Take the following precautions with probe cables:

- Keep free from wheels.
- Do not bend the cable acutely.
- Avoid crossing cables between probes.

Probe orientation

Each probe is provided with an orientation marking. This mark is used to identify the end of the probe corresponding to the side of the image having the orientation mark on the display.

Figure 5-1 Orientation Marking on Probe (Example)



1. Orientation Mark

Labeling

Each probe is labeled with the following information:


- Seller's name and manufacturer
- GE part number
- Probe serial number
- Month and year of manufacture
- Probe designation-provided on the probe grip and the top of the connector housing, so it is easily read when mounted on the system and is also automatically displayed on the screen when the probe is selected.

E7C8L-RS Probe Composition

1. Probe body, including:
 - Probe body
 - Cable
 - Connector
2. Probe Care Card
3. Carrying Case

Illustration for E7C8L-RS probe

Figure 5-2 E7C8L-RS Illustration



1. Knob
2. Twist lock
3. Sterile/Sanitary Outer Sheath


L4-20t-RS Probe Configuration

1. Press **Probe buttons** from **Utility > System**.

Probes and Biopsy

- For each button on the probe handle, select the operation to perform from the drop-down menu.

Figure 5-3 L4-20t-RS Probe Configuration



- Press **Save** to save the configuration.

NOTE

Simultaneously press the button (1, 2) or (3, 4) can enable or disable the probe button functions.

Probe Usage

For details on connecting, activating, deactivating, disconnecting, transporting and storing the probes, refer to *Probes* on page 50 for more information.

Care and Maintenance

Planned maintenance program

The following maintenance schedule is suggested for the system and probes to ensure optimum operation and safety.

Table 5-1 Planned Maintenance Program

Do the Following	Daily	After Each Use	As Necessary
Inspect the Probes	X	X	X
Clean the Probes		X	X
Disinfect Probes		X	X

**CAUTION**

Improper handling can lead to early probe failure and electric shock hazards.

DO follow the specific cleaning and disinfection procedures provided in this chapter and the chemical manufacturers instructions.

Failure to do so will void probe warranty.

**CAUTION**

Transesophageal, endocavitary and intraoperative probes require special handling. Refer to the user documentation enclosed with these probes.

It is recommended to keep a maintenance log and note all probe malfunctions. Follow the maintenance schedule below to ensure optimum operation and safety:

After each use

- Inspect the probe.
- Clean the probe.
- Disinfect the probe.

Before each use

- Inspect the probe.

Environmental Requirements

Probes should be operated, stored, or transported within the parameters outlined below.

For probe environment requirements. See *Environmental Requirements* on page 17 for more information.

**CAUTION**

Ensure that the probe face temperature does not exceed the normal operation temperature range.

Inspecting the probe**CAUTION**

If any damage is found, DO NOT use the probe until it has been inspected and released for further use by a GE HealthCare service representative.

Before each use

1. Inspect the probe.


**CAUTION**

Inspect probes for sharp edge or rough surfaces that could injure sensitive tissue.

Inspect the probe's lens, cable, casing, and connector for cracks, cuts, tears, and other signs of physical damage.

2. Look for damage that might allow liquid into the probe.
3. Test the functionality of the probe.

Table 5-2 Probe parts



After each use

1. Inspect the probe's lens, cable, casing, and connector for cracks, cuts, tears, and other signs of physical damage. Look for any damage that would allow liquid to enter the probe. If any damage is found, do not use the probe until it has been inspected and repaired/replaced by a GE HealthCare Service Representative.
2. Look for damage that might allow liquid into the probe.

Probe Care Cards

The Probe Care Card, supplied with every probe, contains a list of additional processing products that are compatible with probe materials, but GE HealthCare has not validated their efficacy. When using these products, refer to their manufacturer's instructions for use.

The reprocessing instructions provided in this document have been validated with the chemicals specified in *Table 5-5 Chemicals used for Efficacy Validation* on page 305.

A complete list of compatible products is also available on the Ultrasound Probes website. This list updates periodically with new products and probes.

Table 5-3 Website Links

Support Documentation Library Website:
https://www.gehealthcare.com/documentation
Ultrasound Probes Website
https://www.gehealthcare.com/transducers
GE HealthCare chemicals compatibility
https://cleaning.gehealthcare.com/


With the exception of the chemicals listed in *Chemicals Used for Efficacy Validation* on page 304, the chemicals listed on the website have been validated for compatibility only. Therefore, only those agents listed are recommended due to validated efficacy, and any reprocessing performed using chemicals not listed on the website must be validated by the user.

Adequate cleaning and disinfection between patient cases are necessary to prevent disease transmission. All probes must be thoroughly cleaned prior to disinfection. The level of disinfection required is based on patient contact.

- Probes that contact mucosal or non-intact skin require cleaning followed by high-level disinfection by either soaking or use of a trophon EPR or trophon 2.
- Probes that contact intact skin require cleaning followed by intermediate-level disinfection (wipe or spray).
- Verify probe compatibility using *Chemicals Used for Efficacy Validation* on page 304 and the Ultrasound Probe Website listed above.

Table 5-4 Description of Pictogram on Care card

Pictogram	Description
	"ATTENTION" - Consult accompanying documents is intended to alert the user to refer to the operator manual or other instructions when complete information cannot be provided on the label.
	"CAUTION" - Dangerous voltage (the lightning flash with arrowhead) is used to indicate electric shock hazards.
	Biohazard - Patient/user infection due to contaminated equipment. Usage • Cleaning and care instructions • Sheath and glove guidelines
	Symbol that indicates a useful information.
	Pictogram on Probe Care Card: Describes precautions necessary to prevent the risk of disease transmission or infections.
	Ultrasound probes are highly sensitive medical instruments that can easily be damaged by improper handling. Use care when handling and protect from damage when not in use.
	Do not immerse the probe into any liquid beyond the level specified for that probe. Refer to the user manual of the ultrasound system.
	Since there is a possibility of having negative effects on the probe, observe the specified immersing time by the germicide manufacturer strictly. Do not immerse the probe in liquid chemical germicides more than the time prescribed in the care card.
	"Consult accompany document" - Refer to the ultrasound system user manual for important probe care and cleaning instruction.

Pictogram	Description
	This is to illustrate compatible ultrasound coupling gels.
	This is to illustrate compatible cleaners or disinfectants available in spray format (to be used according to instruction from the manufacturers of these products).
	This is to illustrate compatible cleaners or disinfectants available in wipes format (to be used according to instruction from the manufacturers of these products).
	This is to illustrate compatible cleaners or disinfectants available in powder (to be used according to instruction from the manufacturers of these products).
	This is to illustrate compatible cleaners or disinfectants available in liquid format (to be used according to instruction from the manufacturers of these products).
	This is to illustrate compatible automated reprocessors (to be used according to instruction from the manufacturers of these products).



CAUTION

CREUTZFELD-JACOB DISEASE

This device is not indicated for neurological use. Neurological contact on patients with this disease must be avoided. If a device/probe becomes contaminated, there is no adequate means to disinfect it. In this case, the contaminated device/probe must be discarded in accordance with local biologic waste hazard procedures.



CAUTION

Please refer to the probe care card for GE HealthCare approved probe disinfectants.

Chemicals Used for Efficacy Validation

The table below lists the products and intended use (clean, ILD, HLD) that were validated.

Table 5-5 Chemicals used for Efficacy Validation

Product Type	Trade Name	Manufacturer	Minimum Contact Time	Active Ingredient	Validated for Probe
Cleaning (Wipe)	Oxivir® Tb	Diversey	N/A	Hydrogen Peroxide	3Sc-RS, 4C-RS, E8C-RS, 8C-RS, E8Cs-RS, C1-5-RS, IC9-RS, L3-12-RS, L6-12-RS, 9L-RS, 12L-RS, L8-18i-RS, L4-20t-RS, RAB2-6-RS, RIC5-9A-RS, BE9CS-RS
Cleaning (Soak)	Enzol®	Advanced Sterilization Products® (J&J)	1-Minute Soak	Proteolytic Enzymes	3Sc-RS, 6S-RS, 12S-RS, 4C-RS, E8C-RS, 8C-RS, E8Cs-RS, C1-5-RS, IC9-RS, L3-12-RS, L6-12-RS, 9L-RS, 12L-RS, L8-18i-RS, L4-20t-RS, RAB2-6-RS, RIC5-9A-RS, BE9CS-RS, E7C8L-RS
	(Cidezyme®)				
	MetriZyme™	Metrex™			
	Prolystica® 2X Concentrate Presoak & Cleaner	Steris			
Intermediate-level Disinfectant (wipe or liquid)	Oxivir® Tb	Diversey	10-Minute Exposure	Hydrogen Peroxide	3Sc-RS, 4C-RS, 8C-RS, C1-5-RS, L3-12-RS, L6-12-RS, 9L-RS, 12L-RS, L4-20t-RS, RAB2-6-RS
Intermediate-level Disinfectant (wipe)*	Sani-Cloth AF3	PDI	3-Minute Exposure	Quaternary ammonium	12S-RS, 6S-RS
High-level Disinfectant (Soak)	Cidex® OPA (FDA cleared)	Advanced Sterilization Products (J&J)	12-Minute Soak	Ortho-phthalaldehyde	3Sc-RS, 6S-RS, 12S-RS, 4C-RS, E8C-RS, 8C-RS, E8Cs-RS, C1-5-RS, IC9-RS, L3-12-RS, L6-12-RS, 9L-RS, 12L-RS, L8-18i-RS, LK760-RS, RAB2-6-RS, RIC5-9A-RS, BE9CS-RS, E7C8L-RS.
	McKesson OPA/28	McKesson	10-minute soak		

Product Type	Trade Name	Manufacturer	Minimum Contact Time	Active Ingredient	Validated for Probe
High-level Disinfectant (wipe)**	Tristel Trio	Tristel	1-minute soak	Chlorine dioxide	RIC5-9A-RS
	Tristel Duo ULT				

Probe Pre-Cleaning Instructions (Required for All Probes)


The pre-cleaning step is for removal of gel and gross contamination.

NOTE

Manual cleaning is required to ensure the probes are cleaned to the extent necessary for further processing.

1. After each use, disconnect the probe from the ultrasound console. Remove protective sheath from the probe and gently remove all coupling gel from the probe by wiping with soft, lint-free cloth.
2. Wipe the probe with one of the wipes listed in Probe Care Card from the lens, past the strain relief and approximately fifty (50) centimeters down the cable. Wipe the cable with a lint-free cloth dampened with potable water to remove chemical residue. Dispose of the cloth, wipe and gloves in the clinical trash.

Figure 5-4 Cleaning the Probe Cable



NOTE

Use of wipes listed in the Probe Care Card may result in discoloration of the cable.

**WARNING**

Use caution when cleaning the connector. This cable connector should only be cleaned with a slightly dampened cloth or wipe. Exposure to excessive moisture will result in damage to the probe and possibly the ultrasound console. DO NOT wet the connector/console interface surface or labels.

3. After each use, inspect the lens, cable, and housing of the probe. Look for any damage that would allow liquid to enter the probe. If the probe is damaged, do not place it into any liquid (e.g. for disinfection) and do not use it until it has been inspected and repaired/replaced by a GE HealthCare Service Representative.

Probe Manual Cleaning Instructions (Required for All Probes)

1. Put on a clean pair of gloves and fill a sink or basin with warm potable water (30 - 40°C) to a level allowing immersion of the probe up to the immersion line.
2. Prepare the cleaning solution in accordance with the detergent manufacturer's instructions.
3. Immerse the probe in the prepared cleaning solution up to the immersion line and ensure no air bubbles are trapped on the surface. Do not submerge probe beyond the immersion line.

NOTE

Manual cleaning efficacy has been validated using ENZOL® Enzymatic Detergent.

4. Brushing with a clean, soft, nylon bristle brush from the base of the cable strain relief to the distal tip is critical to ensure cleaning and disinfection efficacy.


Figure 5-5 Cleaning the probe using a brush



**CAUTION**

Do not use the brush on the probe lens.

Figure 5-6 Probe Lens Examples



5. Continue brushing the probe for not less than the minimum contact time listed on the detergent manufacturer's label.
6. Visually inspect the probe for soil. Repeat Steps 3 through 5 until all visible soil has been removed from the surface of the probe.
7. Rinse the probe under running warm potable water (30 - 40°C) for not less than 2 minutes. Scrub the surface of the probe with a clean, soft, nylon bristle brush from the base of the cable strain relief to the distal tip.

**WARNING**

User should pay more attention to the protection of probe connector during the cleaning.

**CAUTION**

Take extra care when handling the lens face of the ultrasound transducer. The lens face is especially sensitive and can easily be damaged by rough handling. NEVER use excessive force or abrasive paper products when cleaning the lens face.

8. Visually inspect the device in a well-lit area to ensure all surfaces are free from residual cleaning solution. Repeat Step 7 if visible cleaning solution is observed.
9. Thoroughly dry the probe using a clean lint-free soft and dry cloth or wipe.

Cable and Connector Manual Cleaning

**WARNING**

Use caution when cleaning the connector. This cable connector should only be cleaned with a slightly dampened cloth or wipe. Exposure to excessive moisture will result in damage to the probe and possibly the ultrasound console. DO NOT wet the connector/console interface surface or labels.


1. The cable and connector surfaces can be cleaned with the cleaners or wipes listed in the Probe Care Card.

NOTE

Use of wipes listed in the Probe Care Card may result in discoloration of the cable.

2. Wipe the cable with a lint-free cloth dampened with potable water to remove chemical residue.

Figure 5-7 Cleaning the Probe Cable



Probe Intermediate-Level Disinfection - Spray

For Intermediate-Level Disinfection of surface-contacting probes, choose either the spray or wipe method.


NOTE

Probes that contact only intact skin may be disinfected in this manner. All probes that contact non-intact skin or intact mucous membranes (e.g., endocavitary or transesophageal) require High-Level Disinfection.

1. Put on a new pair of gloves and spray enough disinfectant solution to saturate a new disposable lint-free wipe or cloth.
2. Holding the probe near the strain relief, apply the dampened cloth to the patient contacting lens. Wipe the probe from the lens to the strain relief, slightly rotating the probe after each wiping pass.
3. After the probe has been completely wiped, dampen a second wipe with disinfectant and starting at the probe lens begin wiping the probe in a rotating motion moving down

towards the strain relief. Spray disinfectant directly on the recessed areas and ridges to saturate.

Figure 5-8 Disinfecting the Probe Moving from Lens to Strain Relief



4. Once the probe has been completely wiped, dampen a third wipe with disinfectant and continue wiping the probe as needed to ensure the surface remains wet for the required exposure time. Use as many wipes as needed and respray disinfectant on recessed areas and ridges, to ensure all surfaces remain wet for the minimum required contact time listed in the disinfectant manufacturer's instructions for use.
5. Dry all surfaces of the probe using a soft lint-free wipe or cloth.
6. Saturate a soft lint-free wipe with de-ionized or soft water (remove excess water, wipe should be damp but not dripping) and thoroughly wipe all surfaces of the probe to remove chemical residue. Discard the wipe.
7. A total of three (3) rinses are required. Repeat Step 6 two additional times using new wipes and water.



WARNING

Failure to properly rinse probes with water following disinfection may cause skin irritation.

8. Thoroughly dry all surfaces of the probe using a soft lint-free wipe or cloth, changing wipes/cloths when necessary to ensure the probe is completely dry. Visually inspect the probe to ensure all surfaces are dry. Repeat drying steps if any moisture is visible.
9. After each use, inspect the lens, cable, and housing of the probe. Look for any damage that would allow liquid to enter the probe. If the probe is damaged, do not place it into any liquid (e.g. for disinfection) and do not use it until it has been inspected and repaired/replaced by a GE HealthCare Service Representative.

Probe Intermediate-Level Disinfection - Wipe


1. Put on a new pair of gloves. Holding the probe near the strain relief, apply the wipe to the patient contacting lens. Wipe the probe from the lens to the strain relief, slightly rotating the probe after each wiping pass.
2. After the probe has been completely wiped, use a second wipe and starting at the probe lens begin wiping the probe in a rotating motion moving down towards the strain relief.

Wring the wipe above recessed areas, seams, and ridges to drip disinfectant directly onto these less accessible surfaces.

NOTE


Intermediate-level disinfection efficacy has been validated using Oxivir® Tb Wipes.

Figure 5-9 Disinfecting the Probe Moving from Lens to Strain Relief



3. Once the probe has been completely wiped, use a third wipe and continue wiping the probe as needed to ensure the surface remains wet for the required exposure time. Use as many wipes as needed and drip additional disinfectant on recessed areas and ridges, to ensure all surfaces remain wet for the minimum required contact time listed in the disinfectant manufacturer's instructions for use.

Figure 5-10 Disinfecting the Probe



4. Dry all surfaces of the probe using a soft lint-free wipe or cloth.
5. Saturate a soft lint-free wipe with de-ionized or soft water (remove excess water, wipe should be damp but not dripping) and thoroughly wipe all surfaces of the probe to remove chemical residue. Discard the wipe.

6. A total of three (3) rinses are required. Repeat Step 5 two additional times using new wipes and water.

**WARNING**

Failure to properly rinse probes with water following disinfection may cause skin irritation.

7. Thoroughly dry all surfaces of the probe using a soft lint-free wipe or cloth, changing wipes/cloths when necessary to ensure the probe is completely dry. Visually inspect the probe to ensure all surfaces are dry. Repeat drying steps if any moisture is visible.
8. After each use, inspect the lens, cable, and housing of the probe. Look for any damage that would allow liquid to enter the probe. If the probe is damaged, do not place it into any liquid (e.g. for disinfection) and do not use it until it has been inspected and repaired/replaced by a GE HealthCare Service Representative.

Probe High Level Disinfection - Soak

1. Put on a clean pair of gloves and fill a sink or basin with high level disinfectant diluted in accordance with the disinfectant manufacturers instructions to a level allowing immersion of the probe up to immersion line shown in the user manual and in the illustration below.

NOTE

High-level disinfection efficacy has been validated using Cidex® OPA.

NOTE

Cleaning and disinfection instructions for Transesophageal probes are documented in the Transesophageal Probe Care Card and Probe User Manual.

NOTE

All semi-critical probes that contact non-intact skin or intact mucous membranes require High-Level Disinfection.

NOTE

Handles of semi-critical probes that are not submerged during High-Level Disinfection require low or Intermediate-Level Disinfection to avoid cross contamination.

Figure 5-11 Probe Immersion Levels




Table 5-6 Description of Pictogram for Probe Immersion Levels

Pictogram	Description
	Fluid Level

**WARNING**

Use only high-level disinfectants that are listed in the Probe Care Card enclosed with the probe. In addition, refer to the local / national regulations. Do not steam autoclave or subject the probe to Ethylene Oxide (ETO).

2. Immerse probe in the disinfectant up to the immersion line and ensure no air bubbles are trapped. Ensure the probe remains in the disinfectant for at least the minimum contact time listed in the disinfectant manufacturer's instructions for use.

**WARNING**

Do not immerse the probe in liquid beyond the level specified for that probe.

Never immerse the probe connector or probe adapters in liquid.

The probe should not be exposed to the disinfectant longer than specified to achieve the desired effect.

DO NOT soak or saturate probes with solutions containing alcohol, bleach, ammonium chloride compounds or hydrogen peroxide.

**CAUTION**

Ensure that the probe is suspended. The probe face should not be resting against the tank/basin surface and should be in full contact with the liquid. Carefully place the probe in the basin, taking care not to damage the transducer lens.

Figure 5-12 Probe suspended in disinfectant basin**NOTE**

When performing high-level disinfection of GE HealthCare ultrasound probes with the Trophon® EPR, it is not necessary to disconnect the probe from the ultrasound system. The probe must be inactive (not selected) during the disinfection cycle.

3. Thoroughly rinse the probe by immersing it in a large volume of critical (purified) water for a minimum of 1 (one) minute. Remove the probe and discard the rinse water.

Do not reuse the water. Always use fresh volumes of water for each rinse. Repeat Step 3 two additional times, for a total of 3 (three) rinses.

**WARNING**

Failure to properly rinse probes with water following disinfection may cause skin irritation.

4. Thoroughly dry all surfaces of the probe using a soft, lint-free wipe or cloth, changing wipes' cloths when necessary to ensure the probe is completely dry. Visually inspect the probe to ensure all surfaces are clean and dry. Repeat drying steps if any moisture is visible.
5. After each use, inspect the lens, cable, and housing of the probe. Look for any damage that would allow liquid to enter the probe. If the probe is damaged, do not place it into any liquid (e.g. for disinfection) and do not use it until it has been inspected and repaired/replaced by a GE HealthCare Service Representative.
6. If the probe is not immediately reused, store the probe in a manner that will protect and keep the probe from being recontaminated. This may be accomplished by placing the probe in a storage cabinet with filtered air flow and/or by using a disposable storage cover placed over the probe.

The instructions provided above have been validated to properly prepare GE HealthCare Ultrasound probes for re-use. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

**WARNING****CREUTZFELDT-JACOB DISEASE**

Neurological use on patients with this disease must be avoided. If a probe becomes contaminated, there is no adequate disinfecting means.

**CAUTION**

Ensure that you follow the probe disinfection procedure provided by GE HealthCare.

**CAUTION**

Review the probe care card that is packed with each probe.

**CAUTION**

Please refer to the probe care card for GE HealthCare approved probe disinfectants.

Covering the Transducer using a Sterile, Protective Sheath

**CAUTION**

Protective barriers may be required to minimize disease transmission. Probe sheaths are available for use with all clinical situations where infection is a concern. Use of legally marketed, sterile probe sheaths is mandatory for intra-cavitary and intra-operative procedures.

1. Place an appropriate amount of gel inside the protective sheath and/or on the transducer face.

Probes and Biopsy

2. Insert transducer into sheath, making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing the sheath.

Figure 5-13 Applying the Sheath



- a. Secure the Sheath with a rubber band.
- b. The probe sheath should extend past the end of the probe to the probe's cable.

NOTE

No gel was applied to the probe in this photo.

3. Secure the sheath in place.

NOTE

Failure to use a sheath that fully covers the transducer to the cable strain relief may lead to cross-contamination of the transducer.

4. Inspect the sheath to ensure there are no holes or tears. If the sheath becomes compromised, stop the procedure and replace immediately.

E7C8L-RS Operation

NOTE

Inspect probe for damage or cracks prior to each use. A damaged probe may cause an electrically hazardous condition when coupled to the human body. DO NOT use any probe with damaged lens to scan a patient.

Further more, a damaged probe will not produce a desirable image.

Inspecting before each use or the operation

1. Make a visual inspection of the probe before each use for damage or degradation to the housing, strain relief, lens and seal.
2. Check all surfaces of the probe, especially the patient contact area, to ensure that there are no rough or sharp edges that could damage the sheath or patient tissue.

Tilting the grip

Figure 5-14 Grip adjustment



1. Loosen the screw by turning the knob on the probe body counterclockwise.
2. Set the grip to the desired angle and securely tighten the screw by turning the knob clockwise.



CAUTION

If there is an attempt to tilt the grip or change the grip angle after the sterile sheath is placed on the probe, the sterile sheath may break if it gets caught in the grip. To avoid this, be sure to tilt the grip before putting on the sterile sheath.

System Preparation

Connect the two E7C8L-RS probe connectors to the Ultrasound system's probe ports. Before inserting the connector into the probe port, inspect the probe connector pin. If the pin is bent, do not use the probe until it has been inspected and repaired/replaced by a GE HealthCare Service Representative.

NOTE

Sterile/sanitary sheaths are to be used on the probe during its actual use with patients. Wearing gloves protects the patient and operator.

Patient Preparation for Transrectal Imaging

1. Prepare the patient. An enema is recommended one hour before the exam.
2. Transaxial imaging is best performed with the patient in the supine position.

Change Scan Method

Change scan methods between convex and linear as follows:

- Press the **Probe** key on the touch panel.
- Select probe that corresponds to the probe port used.

The split-screen 2D mode display format may be useful as well.

Patient Scan

1. Scan the patient. The probe handle orientation mark indicates the image scan plane.

2. If necessary, when the probe is in position, attach it to the mechanical stepping device.

Steps for attaching probe to stepper:

- a. Place the probe into clamp.
- b. Tighten the probe clamp.
- c. Attach the needle placement guide.
- d. Align the probe to the needle placement guide.



CAUTION

Ensure that the arrow at the bottom center of the Needle Placement Guide accurately aligns to the center line marked on the probe shaft.

Failure to do so will cause the needle placement not to match the ultrasound display.

3. If a needle placement is to be performed, display the needle placement grid on the image screen.



CAUTION

Before needle insertion, scan the patient to determine the correct puncture depth and site.

Only the sterile/sanitary sheath, rubber band or twist lock and finger cot with rubber band are on the E7C8L-RS probe during the pre-needle placement scanning.


4. When the examination is over:
 - a. Carefully remove the probe.
 - b. Remove the twist lock. Remove and discard the sheath.

- c. Thoroughly clean and disinfect the probe and equipment. Return the probe to its box.

Figure 5-15 Scan Plane Orientation Transaxial View



Figure 5-16 Transrectal Image Orientation Axial View



Probe Safety

Handling precautions



WARNING

Ultrasound probes are highly sensitive medical instruments that can easily be damaged by improper handling. Use care when handling and protect from damage when not in use. DO NOT use a damaged or defective probe. Failure to follow these precautions can result in serious injury and equipment damage.

Electrical shock hazard



ELECTRICAL HAZARD

The probe is driven with electrical energy that can injure the patient or user if live internal parts are contacted by conductive solution:

- **DO NOT** immerse the probe into any liquid beyond the level indicated by the immersion level diagram. Refer to the immersion illustration in the Probe Cleaning Process section. Never immerse the probe connector into any liquid.
- **DO NOT** drop the probes or subject them to other types of mechanical shock or impact. Degraded performance or damage such as cracks or chips in the housing may result.
- Prior to each use, visually inspect the probe lens and case area for cracks, cuts, tears, and other signs of physical damage. **DO NOT** use a probe which appears to be damaged until you verify functional and safe performance. You must perform a more thorough inspection, including the cable, strain relief, and connector, each time you clean the probe.
- Before inserting the connector into the probe port, inspect the probe connector pins. If a pin is bent, do not use the probe until it has been inspected and repaired/replaced by a GE HealthCare Service Representative.
- **DO NOT** kink, tightly coil, or apply excessive force on the probe cable. Insulation failure may result.
- Electrical leakage checks should be performed on a routine basis by GE HealthCare Service or qualified hospital personnel. Refer to the service manual for leakage check procedures.

Mechanical hazards



WARNING

Allergic reactions possible with any probe sheaths and needle used.



WARNING

There is a potential risk that ultrasound imaging may be lost during a critical procedure, such as biopsy scenario. If ultrasound imaging is lost, please restart the device and follow hospital's procedure to take care of patients.



WARNING

Image noise (such as white lines) could be seen during scanning. These can be caused by a short-term external interference source. For example, strong electromagnetic interference equipment (e.g. high frequency electric knife) or larger fluctuations of the power load of other equipment which is connected to the same power grid. Please take mitigation measures, such as relocating or re-orienting the equipment; if the noise is still there, please check and remove the interference source.

**CAUTION**

A defective probe or excessive force can cause patient injury or probe damage:

- Observe depth markings and do not apply excessive force when inserting or manipulating intracavitary probes.
- Inspect probes for sharp edges or rough surfaces that could injure sensitive tissue.
- DO NOT apply excessive force to the probe connector when inserting into the probe port. The pin of a probe connector may bend.

6Tc-RS Probe Thermal Safety

Maintaining a safe thermal environment for the patient has been a design priority at GE. It is generally agreed that in order to avoid damage to body tissues, for long term exposures, tissue contact probe tip temperatures should be less than 42.7°C. The ultrasound scanner incorporates an elaborate thermal safety system which informs the physician of the operating temperature of the probe, and prevents the operative temperature from exceeding given limits. Whenever the 6Tc-RS Probe is plugged into the system and selected, the probe tip temperature is displayed on the system monitor. When the probe tip temperature is over 41.0°C, the temperature displayed on the system monitor turns to red.

NOTE

When the probe tip temperature is higher than 41.0°C, there will be a warning message displays on the screen. User can press **OK** to continue scanning or press **Cancel** to freeze the image.

**WARNING**

The transesophageal probe should be separated from the patient before using the defibrillator.

If the temperature sensor is not working properly when you plug the probe into the system, the probe will not be accepted and scanning will not be possible.

If the probe temperature is over 42.7°C (including probe temperature measurement error), the system will shut down the probe transmitting.

The absolute tolerance of 6Tc-RS probe surface temperature measurement is less than 0.5°C.

High temperature protection levels

The temperature is always displayed on the system monitor. The system has two levels of upper thermal limit: the first high limit is set at 41.0°C, and the second high limit is set at 42.7°C. If the temperature of the probe tip reaches 41.0°C, the temperature display turns red, the system enters freeze mode, and a warning appears on the monitor, asking the user if he/she wishes to continue scanning up to a higher temperature limit. If the response is OK, the scanner will resume scanning. If he/she selects Cancel, or no response is given, the system will remain in Freeze mode. If the temperature reaches 42.7°C, the system will freeze unconditionally. The user will not be allowed to start scanning until the temperature has decreased 0.5°C below the limit where the system entered Freeze mode. To restart scanning, the user must press the Freeze button.

After the probe has been selected, the scan plane positioning system automatically calibrates. This calibration cycle lasts 10 to 15 seconds. After the calibration is complete, the probe temperature sensor is activated, and the probe temperature is displayed.

In case the initialization of the probe fails (no response from the scan plane button after initialization), re-select the probe to repeat the initialization routine.

**CAUTION**

Transesophageal probes require special handling. Refer to the user documentation enclosed with these probes.

Special handling instructions

Using protective sheaths

**CAUTION**

Protective barriers may be required to minimize disease transmission. Probe sheaths are available for use with all clinical situations where infection is a concern. Use of legally marketed, sterile probe sheaths is mandatory for intra-cavitary and intra-operative procedures. Use of legally marketed, sterile, pyrogen free probe sheaths is required for neurological intra-cavity and intra-operative procedures.

Instructions. Custom made sheaths are available for each probe. Each probe sheath kit consists of a flexible sheath used to cover the probe and cable and elastic bands used to secure the sheath.

Sterile probe sheaths are supplied as part of biopsy kits for those probes intended for use in biopsy procedures. In addition to the sheath and elastic bands, there are associated accessories for performing a biopsy procedure which are included in the kit. Refer to the biopsy instructions for the specific probes in the Discussion section of this chapter for further information.

Reordering. To reorder sheaths, please contact your local distributor or the appropriate support resource.

**CAUTION**

Devices containing latex may cause severe allergic reaction in latex sensitive individuals. Refer to FDA's March 29, 1991 Medical Alert on latex products.

**CAUTION**

Do not use pre-lubricated condoms as a sheath. In some cases, they may damage the probe. Lubricants in these condoms may not be compatible with probe construction.

**CAUTION**

DO NOT use an expired probe sheath. Before using probe sheaths, verify whether the term of validity has expired.

Endocavitary Probe Handling Precautions

If the disinfectant solution comes out of the endocavitary probe, please follow the cautions below.

**CAUTION**

Sterile/sanitary sheaths are to be used on the probe during its actual use with patients. Wearing gloves protects the patient and operator.

**CAUTION**

Disinfectant Exposure to Patient: Contact with a disinfectant to the patient's skin or mucous membrane may cause an inflammation. If this happens, refer to the disinfectant's instruction manual.

Disinfectant Exposure from Probe Handle to Patient: DO NOT allow the disinfectant to contact the patient. Only immerse the probe to its specified level. Ensure that no solution has entered the probe's handle before scanning the patient. If disinfectant comes into contact with the patient, refer to the disinfectant's instruction manual.

Disinfectant Exposure from Probe Connector to Patient: DO NOT allow the disinfectant to contact the patient. Only immerse the probe to its specified level. Ensure that no solution has entered the probe's connector before scanning the patient. If disinfectant comes into contact with the patient, refer to the disinfectant's instruction manual.

Endocavitary Probe Point of Contact: Refer to the disinfectant's instruction manual.

Probe handling and infection control

This information is intended to increase user awareness of the risks of disease transmission associated with using this equipment and provide guidance in making decisions directly affecting the safety of the patient as well as the equipment user.

Diagnostic ultrasound systems utilize ultrasound energy that must be coupled to the patient by direct physical contact. Depending on the type of examination, this contact occurs with a variety of tissues ranging from intact skin in a routine exam to recirculating blood in a surgical procedure. The level of risk of infection varies greatly with the type of contact.

One of the most effective ways to prevent transmission between patients is with single use or disposable devices. However, ultrasound transducers are complex and expensive devices that must be reused between patients. It is very important, therefore, to minimize the risk of disease transmission by using barriers and through proper processing between patients.

**CAUTION**

To reduce the risk of infection, ALWAYS clean and disinfect the probe according to the probe specific instructions including probe compatible chemicals between patients to the level appropriate for the type of examination and use FDA-cleared probe sheaths where appropriate.

**CAUTION**

Adequate cleaning and disinfection are necessary to prevent disease transmission. It is the responsibility of the equipment user to verify and maintain the effectiveness of the infection control procedures in use. Always use sterile, legally marketed probe sheaths for intra-cavitary and intra-operative procedures.

For neurological intra-operative procedures, use of a legally marketed, sterile, pyrogen free probe sheath is REQUIRED. Probes for neuro surgical use must not be sterilized with liquid chemical sterilants because of the possibility of neuro toxic residues remaining on the probe.



CAUTION

In order for liquid chemical disinfectants to be effective, all visible residue must be removed during the cleaning process. Thoroughly clean the probe, as described earlier before attempting disinfection.

You **MUST** disconnect the probe from the ultrasound system prior to cleaning/disinfecting the probe. Failure to do so could damage the system.

DO NOT soak probes in liquid chemical disinfectant for longer than is stated by the disinfectant instructions for use. Extended soaking may cause probe damage and early failure of the enclosure, resulting in possible electric shock hazard.



WARNING

Ultrasound transducers can easily be damaged by improper handling and by contact with certain chemicals. Failure to follow these precautions can result in serious injury and equipment damage.

- Do not immerse the probe into any liquid beyond the level specified for that probe. Never immerse the transducer connector or probe adapters into any liquid.
- Avoid mechanical shock or impact to the transducer and do not apply excessive bending or pulling force to the cable.
- Transducer damage can result from contact with inappropriate coupling or cleaning agents:
 - Do not soak or saturate transducers with solutions containing alcohol, bleach, ammonium chloride compounds or hydrogen peroxide
 - Avoid contact with solutions or coupling gels containing mineral oil or lanolin
 - Avoid temperatures above 50°C
- Inspect the probe prior to use and after use for damage or degeneration to the housing, strain relief, lens and seal. Do not use a damaged or defective probe.

Special Reprocessing Instructions for E7C8L-RS

To properly disinfect the E7C8L-RS probe, the probe handle can be disassembled and reassembled so that the entire probe and handle can be immersed in the germicide solution. Allow the germicide to remain in contact with the fully immersed probe, for high level disinfection, according to the germicide manufacturer's recommended time. To remove all germicide residue, final rinsing should be done following the germicide manufacturer's instructions. Remove excess water by shaking and allow to air dry.

Disassemble the probe handles

1. Turn the probe knob counterclockwise until the pin can be removed.

2. Grasp the probe and slide the left piece (the one with GE Logo) towards the arrow direction to separate the two handle pieces.

NOTE

The indicated direction of the arrows on the handles are only for disassemble process; to reassemble the handles, follow the opposite direction.

Figure 5-17 Disassemble the probe handles

**Reassemble the probe handles**

1. Align the probe with the right handle piece (the one without GE logo).

NOTE

Ensure the probe cables are well positioned to allow the mating piece to be slide in.

2. Place the left handle piece on the right handle piece and slide to lock.

Probes and Biopsy

3. Insert the knob and tighten it clockwise till it is well locked.

Figure 5-18 Reassemble the probe



Instruction for reprocessing the probe

1. Disassemble the probe handles.
2. Dispense the RTU cleaning wipe from the wipe canister to wipe the entire probe, handle piece and locking nut. Wrap a wipe around a soft nylon bristle brush to clean the positioning teeth on the handle and probe body.
3. Use a round, soft bristle brush wrapped with a cleaning wipe to clean the through hole on the probe and handle.
4. Ensure the probe remains exposed to the cleaning solution for at least the minimum contact time listed on the enzymatic cleaner label.
5. Rinse under running water for at least 2 minutes, or until no visible signs of cleaner residue.
6. Visually inspect the device in a well-lit area to ensure that no cleaning residue remains on all surfaces. Repeat Step 5 if visible cleaning solution is observed.
7. Use a clean, lint-free soft and dry cloth or wipe to dry the probe thoroughly. Do not use abrasive paper products.
8. Visually inspect the probe for residue, if necessary, repeat steps 2 to 5 until the probe is visibly clean.
9. Put on a clean pair of gloves.
10. Prepare a basin with High-level disinfectant per manufacturer's disinfectant instructions.
11. Immerse the probe, handle and knob in the disinfectant. Ensure that the probe is immersed above the cable strain relief device, and ensure that no air bubbles are trapped. Ensure that the probe remains in the disinfectant for at least the minimum contact time listed in the manufacturer's disinfectant instructions.
12. Rinse the probe, handle and knob thoroughly with running water for at least 2 minutes.
13. Reassemble the probe handles.

Coupling gels



WARNING

Do not use unrecommended gels (lubricants). They may damage the probe and void the warranty. Please refer to the probe care card for GE HealthCare approved probe gels.

Applying

In order to assure optimal transmission of energy between the patient and probe, a conductive gel or couplant must be applied liberally to the patient where scanning will be performed.



CAUTION

Do not apply gel to the eyes. If there is gel contact with the eye, flush eye thoroughly with water.

Precautions

Coupling gels should not contain the following ingredients as they are known to cause probe damage:

- Methanol, ethanol, isopropanol, or any other alcohol-based product
- Mineral oil
- Iodine
- Lotions
- Lanolin
- Aloe Vera
- Olive Oil
- Methyl or Ethyl Parabens (para hydroxybenzoic acid)
- Dimethyl Silicone
- Polyether glycol based
- Petroleum



WARNING

User should select the Non-toxic and nonirritating GEL.

Returning/Shipping Probes and Repair Parts

US Department of Transportation and GE policy requires that equipment returned for service MUST be clean and free of blood and other infectious substances.

When you return a probe or part for service (Field Engineer or customer), you need to clean and disinfect the probe or part prior to packing and shipping the equipment.

Ensure that you follow probe cleaning and disinfection instructions provided in this manual.

This ensures that employees in the transportation industry as well as the people who receive the package are protected from any risk.

Sterile Ultrasound Procedures

ONLY ultrasound gel that is labeled as sterile, is sterile.

Ensure you always use sterile ultrasound gel for those procedures that require sterile ultrasound gel.

Once a container of sterile ultrasound gel is opened, it is no longer sterile and contamination during subsequent use is possible.

Probe Discussion

The Versana Premier/Versana Premier Lotus supports the following types of probes:

- **Convex Array.** Convex Array probes.
- **Linear Array.** Linear Array probes.
- **Sector Phased Array.** Sector Phased Array probes.
- **Volume Probes.** 4D probe: Convex Array.
- **Bi-plane Probe.** Convex and Linear Array probe and Convex and Convex Array probe.



CAUTION

Probes for transvaginal and transrectal applications require special handling. Transvaginal/transrectal examinations and probe insertions should be performed only by personnel with adequate training. Refer to the user documentation enclosed with these probes.

Application

Table 5-7 Probe Indications for Use

Probe Application	3Sc-RS	6S-RS	12S-RS	4C-RS	E8C-RS	8C-RS	E8Cs-RS	C1-5-RS	IC9-RS	L3-12-RS	L6-12-RS	9L-RS	12L-RS	L8-18i-RS	LK760-RS	L4-20t-RS	6Tc-RS	RAB2-6-RS	RIC5-9A-RS	BE9CS-RS	E7C8L-RS
Abdominal	X			X				X		X	X	X	X			X		X			
Fetal / Obstetrical				X	X		X	X	X	X		X						X	X		
Gynecology				X	X		X	X	X									X			
Vascular/ Peripheral Vascular	X	X	X	X				X		X	X	X	X	X		X					
Urology				X	X		X	X	X									X	X	X	
Pediatric [1]		X	X	X		X		X		X	X	X	X								
Small Parts (includes breast, testes, thyroid)										X	X	X	X	X		X					
Musculoskeletal Conventional [2]				X		X		X		X	X	X	X		X	X					
Musculoskeletal Superficial										X	X	X	X	X		X					
Cardiac Adult	X																X				
Cardiac Pediatric	X	X	X			X															

Probe Application	3Sc-RS	6S-RS	12S-RS	4C-RS	E8C-RS	8C-RS	E8Cs-RS	C1-5-RS	IC9-RS	L3-12-RS	L6-12-RS	9L-RS	12L-RS	L8-18i-RS	LK760-RS	L4-20t-RS	6Tc-RS	RAB2-6-RS	RIC5-9A-RS	BE9CS-RS	E7C8L-RS
Thoracic/Pleural	X			X				X		X	X	X	X			X					
Transcranial	X	X	X			X															
Transrectal					X		X		X										X	X	X
Transvaginal					X		X		X										X		
Transesophageal																	X				
Interventional Guidance																					
Tissue biopsy	X			X	X		X	X	X	X	X	X	X			X	X	X	X	X	X
Fluid drainage	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Vascular										X	X	X	X	X	X	X					
Non-Vascular access	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
[1]: Pediatric application is applicable for Neonatal, Infants, Children and Adolescents, and the age limit ranges from birth up to 21. (Reference: 'Guidance for Industry and FDA Staff: Premarket Assessment of Pediatric Medical Devices' dated March 24, 2014.) [2]: Musculoskeletal Conventional is equal to Musculoskeletal General.																					

Features

Table 5-8 Probe modes of operation and Features

Probe modes of operation and Features	C1-5-RS	4C-RS	8C-RS	E8C-RS	E8Cs-RS	BE9CS-RS	E7C8L-RS	L3-12-RS	L6-12-RS	9L-RS	12L-RS	L8-18i-RS	LK760-RS	3Sc-RS	6S-RS	12S-RS	RAB2-6-RS	RIC5-9A-RS	L4-20t-RS	IC9-RS	6Tc-RS
B Mode	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
M Mode	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
PW Mode	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
CW Mode														X	X	X					X
CF Mode	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Color M Mode														X	X	X					X
PDI	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X

Probe modes of operation and Features	C1-5-RS	4C-RS	8C-RS	E8C-RS	E8Cs-RS	BE9CS-RS	E7C8L-RS	L3-12-RS	L6-12-RS	9L-RS	12L-RS	L8-18i-RS	LK760-RS	3Sc-RS	6S-RS	12S-RS	RAB2-6-RS	RIC5-9A-RS	L4-20t-RS	IC9-RS	6Tc-RS
Directional PDI	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Combined Modes	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Harmonic Imaging	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
TVI														X	X	X					X
TVM														X	X	X					X
TVD														X	X	X					X
Static 3D/4D																	X	X			
BFlow	X	X						X	X	X	X								X		
BFC	X	X						X	X	X	X								X		
Contrast	X	X						X												X	
Elastography	X	X			X			X	X	X	X									X	
Anatomical M	X	X	X											X	X	X	X				X
Curved Anatomical M														X	X	X					X
Easy 3D	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X			X	X	
Advanced 3D	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X			X	X	
TUI																	X	X			
TVI w/ QAnalysis														X	X	X					X
Logiq View	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
CrossXBeam	X	X	X	X	X	X	X	X	X	X	X	X	X				X	X	X	X	
SRI HD	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Auto IMT								X	X	X	X	X							X		
Stress Echo														X							
Auto EF														X							
Virtual Convex								X	X	X	X	X	X	X	X	X			X		
Sonobiometry	X	X															X				
Vocal																	X	X			

Probes and Biopsy

Probe modes of operation and Features	C1-5-RS	4C-RS	8C-RS	E8C-RS	E8Cs-RS	BE9CS-RS	E7C8L-RS	L3-12-RS	L6-12-RS	9L-RS	12L-RS	L8-18i-RS	LK760-RS	3Sc-RS	6S-RS	12S-RS	RAB2-6-RS	RIC5-9A-RS	L4-20t-RS	IC9-RS	6Tc-RS
Breast Productivity and ACR BI-RADS								X	X	X	X								X		
Thyroid Productivity and ACR TI-RADS								X	X	X	X	X							X		
Auto Bladder	X	X		X	X	X	X											X		X	
Needle Recognition	X	X					X	X	X	X	X								X		
Auto TGC	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
ATO	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
ASO	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Whizz	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Whizz Color Flow	X	X						X	X		X										
Whizz Easy Style	X	X						X									X				
Touch LGC	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Digital TGC	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
V-Live 2.0																	X	X			
Whizz RenderLive																	X	X			
Whizz Follicle																	X	X			
V-Zoom	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Strain														X	X	X					X
Strain Rate														X	X	X					X
Whizz Label	X	X																			
Intensity Ratio	X	X	X					X	X	X	X			X			X				
Cine	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
ECG														X	X	X					X
ARC Li-RADS	X	X						X	X	X	X			X			X		X		

Probe modes of operation and Features	C1-5-RS	4C-RS	8C-RS	E8C-RS	E8Cs-RS	BE9CS-RS	E7C8L-RS	L3-12-RS	L6-12-RS	9L-RS	12L-RS	L8-18i-RS	LK760-RS	3Sc-RS	6S-RS	12S-RS	RAB2-6-RS	RIC5-9A-RS	L4-20t-RS	IC9-RS	6Tc-RS
Breast Care								X	X	X	X								X		
Follow-up Tool	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Brachytherapy							X														
Probe Check	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
eDelivery	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Imaging Insights	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Insite AutoConnect	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
DICOM	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Scan Coach	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Scan Assistant	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
My Trainer	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Whizz Report	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Simplified Workflow	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Multi-Touch	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Standby	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Whizz Note	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Digital Keyboard	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X

Specifications

Table 5-9 System Probe Definitions

Probe Designation	Center Image Frequency (MHz)	Frequency Range (MHz)	Doppler Frequency Range (MHz)
3Sc-RS	2.75 MHz $\pm$ 10%	1.7~4.0	1.7~3.3
6S-RS	4.5 MHz $\pm$ 10%	3.0~7.0	3.0~4.5
12S-RS	7.75 MHz $\pm$ 10%	4.2~12.0	4.2~6.7

Probe Designation	Center Image Frequency (MHz)	Frequency Range (MHz)	Doppler Frequency Range (MHz)
4C-RS	3.10 MHz $\pm$ 10%	2.0~5.0	2.0~3.3
E8C-RS	6.5 MHz $\pm$ 20%	4.2~10.0	4.2~6.3
8C-RS	6.5 MHz $\pm$ 20%	4.2~10.0	4.2~6.3
E8Cs-RS	6.5 MHz $\pm$ 20%	4.0~10.0	4.0~6.0
C1-5-RS	3.4 MHz $\pm$ 20%	2.0~5.0	2.0~3.6
IC9-RS	6.25 MHz $\pm$ 0.35	3.6~9.0	3.6~6.7
L3-12-RS	6.575 MHz $\pm$ 1	3.6~12.0	3.6~9.1
L6-12-RS	7.75 MHz $\pm$ 20%	4.0~13.0	4.0~6.0
9L-RS	5.25MHz $\pm$ 20%	4.0~10.0	4.0~5.0
12L-RS	7.5 MHz $\pm$ 20%	4.2~13.0	4.2~7.7
L8-18i-RS	9.5 MHz $\pm$ 20%	6.7~18.0	6.7~10.0
LK760-RS	7.15 MHz $\pm$ 10%	3.5~10.0	3.5~5.0
L4-20t-RS	9.1 MHz $\pm$ 20%	5.0~20.0	5.0~10.0
6Tc-RS	4.8 MHz $\pm$ 0.4	3.0~8.0	3.0~5.0
RAB2-6-RS	3.3 MHz $\pm$ 10%	2.0~6.0	2.0~4.0
RIC5-9A-RS	6.5 MHz $\pm$ 10%	4.2~10.0	4.2~6.3
BE9CS-RS	7.5 MHz $\pm$ 20%	4.0~10.0	4.0~6.3
E7C8L-RS (-L: Linear Transducer)	6.5 MHz $\pm$ 0.5	5.0~10.0	5.0~6.7
E7C8L-RS (-C: Convex Transducer)	6.5 MHz $\pm$ 0.5	5.0~10.0	5.0~6.7

Slice Thickness Specification


Table 5-10 System Probe Definitions

Probe Designation	Probe Slice Thickness Specification [mm]
3Sc-RS	$\leq$ 16.0mm
6S-RS	$\leq$ 8.0mm
12S-RS	$\leq$ 8.0mm
4C-RS	$\leq$ 12.0mm
E8C-RS	$\leq$ 13.0mm
8C-RS	$\leq$ 13.0mm

Probe Designation	Probe Slice Thickness Specification [mm]
E8Cs-RS	$\leq 13.0\text{mm}$
C1-5-RS	$\leq 12.0\text{mm}$
IC9-RS	$\leq 13.0\text{mm}$
L3-12-RS	$\leq 8.0\text{mm}$
L6-12-RS	$\leq 8.0\text{mm}$
9L-RS	$\leq 14.0\text{mm}$
12L-RS	$\leq 8.0\text{mm}$
L8-18i-RS	$\leq 8.0\text{mm}$
LK760-RS	$\leq 12.0\text{mm}$
L4-20t-RS	$\leq 8.0\text{mm}$
6Tc-RS	$\leq 8.0\text{mm}$
RAB2-6-RS	$\leq 12.0\text{mm}$
RIC5-9A-RS	$\leq 13.0\text{mm}$
BE9CS-RS	$\leq 13.0\text{mm}$
E7C8L-RS (-L: Linear Transducer)	$\leq 8.0\text{mm}$
E7C8L-RS (-C: Convex Transducer)	$\leq 13.0\text{mm}$

Probe Illustration

Figure 5-19 Probe illustration - Example



1. Connector: Socket to connect console receptacle for transmitting electrical signals.

Probes and Biopsy

2. Cable: Transmit electrical signals between transducer and connector, wrapped with soft plastics.
3. Housing: A molded hard plastic part to housing transducer and forming the outer shape of a probe.

Sector Probe

Table 5-11 Sector Probe Illustration

Probe	Illustration	Probe	Illustration
3Sc-RS		6S-RS	
12S-RS			

Convex Probe

Table 5-12 Convex Probe Illustration

Probe	Illustration	Probe	Illustration
4C-RS		E8C-RS	
8C-RS		E8Cs-RS	
C1-5-RS		IC9-RS	

Linear Probe

Table 5-13 Linear Probe Illustration

Probe	Illustration	Probe	Illustration
L3-12-RS		L6-12-RS	
9L-RS		12L-RS	
L8-18i-RS		L4-20t-RS	
LK760-RS			

TEE (transesophageal) Probe

Table 5-14 TEE (transesophageal) Probe Illustration

Probe	Illustration	Probe	Illustration
6Tc-RS			

Mechanical 4D Probe

Table 5-15 Mechanical 4D Probe Illustration

Probe	Illustration	Probe	Illustration
RAB2-6-RS		RIC5-9A-RS	

Bi-plane Probe

Table 5-16 Bi-plane Probe Illustration

Probe	Illustration	Probe	Illustration
BE9CS-RS		E7C8L-RS	

Biopsy Special Concerns

Precautions Concerning the Use of Biopsy Procedures

**WARNING**

Only qualified and trained physicians or sonographers could do TEE probe scanning and biopsy procedures.

**WARNING**

Do not freeze the image during a biopsy and TEE procedure. The image must be live to avoid a positioning error.

Biopsy guide zones are intended to assist the user in determining optimal probe placement and approximate the needle path. However, actual needle movement is likely to deviate from the guideline. Always monitor the relative positions of the biopsy needle and the subject mass during the procedure.

**CAUTION**

The use of biopsy devices and accessories that have not been evaluated for use with this equipment may not be compatible and could result in injury.

**CAUTION**

The invasive nature of biopsy procedures requires proper preparation and technique to control infection and disease transmission. Equipment must be cleaned as appropriate for the procedure prior to use.

- Follow the probe cleaning and disinfection procedures and precautions to properly prepare the probe.
- Follow the manufacturer's instructions for the cleaning of biopsy devices and accessories.
- Use protective barriers such as gloves and probe sheaths.
- After use, follow proper procedures for decontamination, cleaning, and waste disposal.

**CAUTION**

Improper cleaning methods and the use of certain cleaning and disinfecting agents can cause damage to the plastic components that will degrade imaging performance or increase the risk of electric shock.

**WARNING**

NEVER reuse the TR5° disposable biopsy guide attachment and Disposable sterile Ultra-Pro II Needle guide kits.

**WARNING**

Please wear sterilized gloves during the biopsy procedure, meanwhile do not operate the console and any non-sterilized area in case any cross infection.

**WARNING**

Please wear sterilized gloves during scanning with Transvaginal & Transrectal probe and the probe shall with one sterilized rubber condoms in case any cross infection.

Preparing for a Biopsy

Displaying the Guide zone

Activate the Biopsy Kit by selecting it from the B-Mode.




WARNING

There is potential risk in biopsy, brachytherapy and TEE exam scenario without imaging guidance.

The available biopsy options appear when Biopsy Kit is selected. There are fixed and adjustable angle biopsy kits and plastic/disposable and reusable biopsy guides available with the ultrasound system depending on the probe. Select the desired biopsy kit.

Rotate Biopsy Kit rotary to turn it off.

Figure 5-20 Biopsy Guide zones for the 4C-RS Probe



The biopsy guide zone represents a path of the needle. The dots which make up the guide zones is the depth readout where:

- Yellow represent 1 cm increments.
- Red represents 5 cm increments.

The display should be carefully monitored during a biopsy for any needle deviation from the center line or guide zone.

The angular accuracy of the biopsy display shall be ± 5 degrees (angle between the center line and the two biopsy guide lines and the needle is inserted).

Before scanning, verify the needle can be visualized within the imaging plane. User appropriate needle length to reach target area.

The Biopsy Guide zone adjusts along with image adjustments, such as image inversion/rotations, zoom and depth changes.

Probes and Biopsy

The needle may vary from the center line or guide zone for various reasons:

- Needle barrel to needle clearance or strength.
- Bracket manufacturing tolerance.
- Needle deflection due to tissue resistance.
- Needle size chosen. Thinner needles may deflect more.

Table 5-17 Biopsy Guide Availability

Probe	Fixed Angle	Multi-Angle		
		MBX1 (cm)	MBX2 (cm)	MBX3 (cm)
4C-RS		4.0	6.0	10.0
E8Cs-RS	TR5: 85 Degree (15.3 cm) RU: 90 Degree			
3Sc-RS		4.2	5.7	8.2
L6-12-RS		1.5	2.5	3.5
12L-RS		1.5	2.5	3.5
RAB2-6-RS		4.0	6.0	8.0
E8C-RS	TR5: 85 Degree (15.3 cm) RU: 90 Degree			
RIC5-9A-RS	RIC59_Single:90 Degree RIC59_RU:90 Degree			
BE9CS-RS	0.6	4.0	6.0	8.0
C1-5-RS		4.0	6.0	8.5
9L-RS		4.0	5.5	7
L3-12-RS		1.5	2.5	3.5
IC9-RS	TR5: 89 Degree			

Table 5-18 Biopsy Guide Availability for E7C8L-RS

Probe	Multi-Angle
E7C8L-RS (-L: Linear transducer)	MBX1, MBX2, MBX3, MBX4, MBX5, MBX6, MBX7, MBX8
E7C8L-RS (-C: Convex transducer)	MBX, GBX1, GBX3

Table 5-19 Biopsy Guide Availability for L4-20t-RS

Probe	Fixed Angle	Multi-Angle				
		MBX1 (degree)	MBX2 (degree)	MBX3 (degree)	MBX4 (degree)	MBX5 (degree)
L4-20t-RS		60	48	38	29	23

**DANGER**

Failure to match the guide zone displayed to the guide may cause the needle to track a path outside the zone.

It is extremely important that when using the adjustable angle biopsy guides, the angle displayed on the screen matches the angle set on the guide, otherwise the needle will not follow the displayed guide zone which could result in repeated biopsies or patient injury.

Preparing the Biopsy Guide Attachment

Convex, Sector and Linear probes have optional biopsy guide attachments for each probe. The guide consists of a non-disposable bracket to attach to the probe, disposable needle clip to attach to the bracket, sheath, gel (sterile gel if necessary) and disposable needle barrels.

The disposable needle barrels are available for a variety of needle sizes.

**CAUTION**

Please refer to the manufacturer's instructions included in the biopsy kit.

The bracket is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure bracket is properly cleaned, sterilized or disinfected before each use.

Disposable components are packaged sterile and are single-use only. Do not use if integrity of packing is violated or if expiration date has passed.

Fixed Needle Biopsy Guide Assembly

**WARNING**

DO NOT use the needle with the catheter (soft tube). There is a possibility of breaking the catheter in the body.

**CAUTION**

Before inserting the needle, scan the patient to determine the correct puncture depth and site. Only the sterile/sanitary sheath and rubber band are on the probe during the pre-needle placement scanning.

Preparation


To prepare the endocavitary probe for use:

1. Remove the probe from the box and carefully examine it for any damage.

Probes and Biopsy

2. If the biopsy guide is to be attached, use the filling removal tool to clean out the attachment area on the probe head.

Figure 5-21 Attachment Filling Removal



- a. Probe Head
- b. Attachment
- c. Filling Removal Tool

3. Clean, then disinfect the probe.

NOTE

Ensure that protective gloves are worn.

Installing the sheath

To install the sheath:

1. Remove the sheath from its package. Do not unroll the sheath.

NOTE

Remember to rinse all sanitary probe sheaths of powder before placing on the probe. Powder can degrade the displayed image.

2. Place an adequate amount of ultrasound gel inside the sheath tip (the gel is between the sheath inner surface and the probe aperture).


NOTE

Ensure that only acoustic coupling gel is used for this purpose.

3. Place the sheath tip over the probe aperture and then pull the sheath end toward the probe handle.

4. Inspect the sheath for nicks, cuts or tears.

Figure 5-22 E8C-RS Probe with Sheath



- a. Probe Handle
- b. Probe Body
- c. Sanitary Sheath

5. Rub a finger over the tip of the probe to ensure all air bubbles have been removed.

Biopsy Guide Preparation

1. If a biopsy is to be performed, snap the metal or plastic biopsy guide on to the probe over the sheath.



CAUTION

Patient injury or repeated biopsies may result. The needle placement will not be as intended if the needle guide is not properly seated and secure.

Figure 5-23 Disposable Biopsy Guide 5 degree Angle





Figure 5-24 Reusable Biopsy Guide



- a. Fix with a screw
2. Place an adequate amount of ultrasound gel on the gel-filled sheath tip's outer surface.
3. Ensure the guide is properly seated and secure by pushing forward on the needle insertion end of the guide until the attachment node is firmly in place in it's hole.

Multi Angle Biopsy Guide Assembly



WARNING

DO NOT attempt to use the biopsy bracket and needle guide until the manufacturer's instructions, provided with the biopsy bracket and needle guide in the kit, have been read and thoroughly understood.

1. Scan the patient and identify the target for biopsy. Move the probe to locate the target to the center of the image. Enable the system biopsy guide zone and try guide zone angles MBX1 to MBX3 to decide the best angle setting for needle path.

Figure 5-25 Example




2. Pull up on the knob (a) to freely move the needle guide attachment. Align the knob with the selected position of the needle guide attachment.
Push the knob down (b) into the desired slot to secure the angle position of the needle guide attachment.

Figure 5-26 Pull up and push down the knob



3. Fit a convex of the biopsy bracket (a) in a concave of the probe (b).


Figure 5-27 Probe/Bracket Alignment



Probes and Biopsy


Hold the side (a) and tuck down the needle guide side (b) until it clicks or locks in place.

Figure 5-28 Probe/Multi-angle Bracket Alignment 2




4. Place an adequate amount of coupling gel on the face of the probe.
5. Place the proper sanitary sheath tightly over the probe and biopsy bracket. Use the rubber bands supplied to hold the sheath in place.

Figure 5-29 Applying Sanitary Sheath



6. Snap the needle guide onto the biopsy guide bracket.

Figure 5-30 Snap the needle guide



7. Push the locking mechanism towards the bracket to secure the lock (a). Make sure the needle guide is firmly attached to the bracket.

Figure 5-31 Lock the Needle guide



8. Choose the desired gauge (size) needle barrel. Twist it back and forth to remove it from the plastic tree.

Figure 5-32 Needle Barrel



9. Place the needle barrel into the needle clip with the desired gauge facing the needle clip and snap into place.

Figure 5-33 Needle Barrel Installation



Remove the biopsy guide


Hold the other side and push out the needle clip attachment side.

Figure 5-34 Remove the biopsy guide



Releasing the needle

According to the following procedure, you remove the needle from a probe and an assembly without moving the needle.



1. Push the knob portion of a sleeve in the direction of the arrow.
2. The needle is released from the assembly.
3. Push the probe and the assembly in the direction of the larger arrow to remove the needle.

4D Biopsy Guide Assembly - Representative Example

Figure 5-35 Mounting the Biopsy Needle Guide to the 4D Probe



1. Fit the biopsy guide bracket onto the probe.
2. Fix the bracket by locking the frame on the opposite side of the needle guide.
3. Snap the needle guide onto the biopsy guide bracket, refer to *Multi Angle Biopsy Guide Assembly* on page 346 for detailed information.

NOTE

Needle guide sterilization with autoclave is recommended if possible for needle guide sterilization.

4D Probe Biopsy Needle Path Selection

To select the needle path and verify that the path of the needle is accurately indicated within the guide zone on the system monitor, perform the following before use:

1. Properly install the bracket and biopsy guide.
2. Scan in a container filled with water (47°C).
3. Select **Biopsy kit**. The available biopsy options from the Touch Panel.

Select the biopsy guide zone where the needle echo passes through the center of the guide zone. Use the selected biopsy guide zone when performing the biopsy

Biopsy Needle Path Verification

To verify that the path of the needle is accurately indicated within the guide zone on the system monitor, perform the following:

- Properly install the bracket and biopsy guide.
- Scan in a container filled with water (47°C).
- Display the biopsy guide zone on the monitor.
- Ensure that the needle echo falls within the guide zone markers.

The Biopsy Procedure

**WARNING**

Biopsy procedures must only be performed on live images.

1. Place coupling gel on the scanning surface of the probe/sheath/biopsy guide assembly.
2. Activate the biopsy guide zone on the system through the B-Mode menu. When using multi-angle guides, ensure that the proper guide zone angle is displayed.
3. Scan to locate the target. Center the target in the electronic guide zone path.

NOTE

Enabling color flow would allow for visualization of the vascular structure around the area to be biopsied.

4. Place the needle in the guide between the needle barrel and needle clip. Direct it into the area of interest for specimen retrieval.

Post Biopsy

When the biopsy is complete, remove the needle barrel, needle clip and probe sheath. Properly dispose of these items in accordance with current facility guidelines.

Clean and disinfect the probe per the instructions provided in this manual.

The biopsy bracket can be cleaned and disinfected in a recommended disinfecting agent and reused.

**CAUTION**

When the biopsy needle guide kit is opened, all parts must be discarded after the procedure whether they have been used or not.

Chapter 6

Advanced Features

In this section

Whizz Easy Style.	355
Whizz CF.	356
B Flow / B Flow Color (Option).	358
Stress Echo (Option).	361
ECG Option.	383
Using 3D (Option).	392
Elastography (Option).	404
Contrast Imaging (Option).	408
Needle Recognition (Option).	438
Breast Care.	443
Breast Productivity Package.	446
Thyroid Productivity Package.	452
Scan Assistant.	460
My Trainer.	468
Follow-up Tool.	470
QAnalysis (Option).	476
Whizz Label (Not applicable for China).	488
Advanced urology procedures.	492
Strain Imaging.	503
Strain rate.	506
Using 4D.	508
Whizz RenderLive.	533
Whizz Follicle.	535
Whizz Note.	540

Advanced Features

VCO for customer. 544

Digital Expert Connect. 558

Intensity Ratio. 559

US LI-RADS (ACR). 564

Whizz Report. 568

Scan Coach. 582

Whizz Easy Style

Whizz Easy Style is a feature intended to provide different image style/presentation for user to choose. It provides different appearance on tissue smoothness, contrast, border delineation and tissue differentiation. The existing GE HealthCare image style is still available for the user to choose.

Whizz Easy Style is only available on specific probes with specific presets. Refer to *Table 5-8 Probe modes of operation and Features* on page 330 for detail information.


Table 6-1 Whizz Easy Style on each probe

Probe	Presets
4C-RS	Abdomen, Renal, OB-2/3, OB-1, Fetus Echo
C1-5-RS	Abdomen, Renal, OB-2/3, OB-1, Fetus Echo
RAB2-6-RS	OB-2/3, OB-1, Fetus Echo
L3-12-RS	Nerve, Shoulder, Knee, MSK GEN, MSK Sup and Thyroid

Activating

Whizz Easy Style is enabled by clicking **Whizz Easy Style** on Touch Panel, and then select **Level 1** or **Level 2**. To exit, press **Whizz Easy Style** on the Touch Panel, and then select **OFF** to disable it.

Figure 6-1 Whizz Easy Style on touch panel



Whizz CF

Whizz CF optimizes the color flow by adjusting color flow gain and frequency automatically when activated by user. Auto gain adjusts the color flow gain according to noise of color image and depth of color ROI. Auto frequency adjusts the frequency with the depth change automatically. This may provide better resolution when ROI is in the near field, or penetration when ROI in the far field. Whizz color flow mode will be turned off automatically if the user adjusts color gain or frequency manually. Adjusting other parameters will not turn this feature off.

NOTE

If user is not satisfied with Whizz color, color gain and frequency can be adjusted or Whizz color can be turned off.


Using Whizz CF

Whizz CF only supports available on specific probes with specific modes. Refer to *Table 5-8 Probe modes of operation and Features* on page 330 for detail information.

To activate Whizz CF,


1. Select **Utility > Whizz > Whizz CF**.

Figure 6-2 Whizz CF configuration



2. Press **Whizz** key under CF mode, the backlight of **Whizz** key turns to green and the Whizz parameter on the scanning screen will display as **On**.

Figure 6-3 Whizz CF On

The image shows a scanning screen with a black background and white text. At the top left, the letters 'CF' are displayed in a larger font. Below this, a list of parameters and their values is shown. The 'Whizz' parameter at the bottom is highlighted with a red rectangular border. The parameters and values are: Frq 2.0, Gn 29.5, L/A 2/3, PRF 1.1, WF 198, S/P 4/10, and Whizz On.

CF	
Frq	2.0
Gn	29.5
L/A	2/3
PRF	1.1
WF	198
S/P	4/10
Whizz	On

To exit **Whizz CF**, press **Whizz** key on the control panel or adjust CF gain or frequency manually.

If user doesn't adjust CF gain/frequency manually when **Whizz CF** is on, the gain and frequency will return to the value before **Whizz CF** is activated. If user adjusts CF gain/frequency manually when **Whizz CF** is on, **Whizz CF** will be turned off and the gain and frequency will stay as the adjusted value.

B Flow / B Flow Color (Option)

B Flow

B-Flow is intended to provide a more intuitive representation of non-quantitative hemodynamics in vascular structures.

B-Flow is digitally-encoded Ultrasound, using digital codes to enhance weak signals from small particulate reflectors (blood flow) and suppress signals from strong reflectors (tissue). Flow and tissue are displayed simultaneously without threshold decision and overlay.

All B-Mode measurements are available with B-Flow active: depth, distance along a straight line, % stenosis, volume, trace, circumference, and enclosed area.

Presetting

Preset the default B Flow Mode, either B Flow or B Flow Color via **Utility > Imaging > BF/BFC** button. Specify the default to be either B Flow or B Flow Color.

Activating

To activate/deactivate B-Flow, press **B-Flow** button on the touch panel or user-defined hard key. Doppler Mode is available while in B-Flow; however, M-Mode and Color Flow/PDI Modes are not available.

Using B-Flow

To optimize the image:

Adjust the frequency, display depth, and focal zone location based on the patient body type and anatomy of interest. Adjust Sensitivity/PRI and Background setting as needed (see below).

Adjust the remaining Imaging parameters and presets as needed; functionally is the same as B-mode when in B-Flow mode.

Scanning Hints

B-Flow provides an intuitive view of blood flow, acute thrombosis, soft plaque, small vessel perfusion, and high grade stenosis. Compared to Color Doppler, it does not display bleeding, blooming, or aliasing artifacts.

The greater the speed, the better the image scatter density and size. If the scan direction is the same as the flow direction, then the image scatter is elongated; if the scan direction is the opposite as the flow direction, then the image scatter is tighter. Therefore, have the scan direction opposite to that of flow direction. Switch the way you hold the probe, with the probe orientation marker inferior to maintain correct orientation on the monitor. Flow starts from where the focal zone is located.

Turn off the background when imaging the kidney, liver, and spleen. Keep the focal zone as close to flow as possible. It is beneficial to narrow the sector width and increase the frame rate.

Benefits

Compared to Color Doppler mode, B-Flow provides better spatial and temporal resolution, displays blood flow in the entire image, i.e. NO ROI, and is not angle dependent as it does not use the Doppler Principle. B-Flow is therefore a more realistic (intuitive) representation of flow information, allowing you to view both high and low velocity flow at the same time.

Affect on other controls

When you activate B-Flow, the system remembers the imaging parameters set while in B-Mode. When you optimize the frame rate via Line Density, you compromise the resolution and when you optimize the resolution, you compromise the frame rate. B-Flow is not available in 3DView.

Bioeffects

Using B-Flow Frequency, Focus Position/Number, Sensitivity/PRI, Line Density, Visualization, and Power Output may change the TI and/or MI. Observe the output display for possible effects.

Flow Type

There are two type of flow type available in B Flow: Low or High.

- Low= Most sensitivity for slow flow (venous, small parts). This setting has the slowest frame rate.
- High= Most sensitivity for fast flow. Provides the highest frame rate and better detection of flow dynamics. Use whenever possible to capitalize on the high frame rate.

Accumulation

Description

Accumulation enhances the flow in an image; ideal to capture dynamic flow in a still picture.

Values

Off - Infinite. Infinite provides the same result as applying CINE Capture to a B-Flow CINE clip.

Benefit

Accumulation detects the maximum signal and holds it (accumulates it) for the level specified (Off - Infinite).

Background

Description

In the Background ON mode, both of B-Flow and B-mode are displayed in one image.

Value

On or Off.

Sensitivity/PRI

Sensitivity/PRI (Pulse Repetition Interval) is proportional to the time interval between the pulses sent to develop the B-Flow image.

In general, a larger value is recommended for slow flow as slow flow detection requires more time separation between pulses so the system can detect the difference in flow profile. However, a larger value could cause bar artifacts on the image. Therefore, it is suggested to not increase the PRI value more than needed. A small value of PRI should be used when the interest is in fast flow only, e.g. viewing a jet in a stenosis case, where the jet is of interest.

NOTE

Sensitivity/PRI is Probe and Model dependent.

B-Flow Color (BFC)

B-Flow Color is intended to provide a B-Flow representation with colorized flow and background B-Mode image.

Although the BFC is based on B-Flow technology, the BFC image is processed and made by the Color Flow processor and therefore has both the advantages of B-Flow and Color Flow. Consequently, imaging parameters and presets are functionally the same as Color Flow/PDI.

Activating

To activate B-Flow Color,

- While in B-Mode, press B-Flow on the touch panel or user-defined hard key to activate B-Flow. Then press Color on the Touch Panel.
- Adjust the parameters to get the best image.

Enhance (B-Flow Color)

Enhance provides a range of choices for B-Flow Color image quality.

Values

Enhance “Off” is the original B-Flow Color setting which gives you the dynamic appearance of the flow with high frame rate.

Enhance “On” adds sensitivity and stability/continuity to the flow’s appearance, along with lower frame rates than Enhance “Off.”

Scale

Same as Scale in Color Flow/PDI Mode.

Values

Increase for higher flow states.

Decrease to display smaller vessels and slower flow.

Stress Echo (Option)

The Versana Premier/Versana Premier Lotus ultrasound system provides an integrated stress echo package, with the ability to perform image acquisition, review, image optimization, and wall segment scoring and reporting for a complete, efficient stress echo examination.

The stress package provides a protocol template for the two types of stress exams (exercise and pharmacological stress).

In addition to preset factory protocol templates, templates can be created or modified to suit your needs.

You can define various quad screen review groups, in any order and combination, that will suit your normal review protocol.

When reviewing stress examination images, the images are viewed at their original image quality, and different post-processing and zoom factors may be applied to the images under review for effective image optimization.

The protocol template may be configured for Continuous Capture.

A stress echo examination consists of three steps:

- Selection of a stress test protocol template
- Image acquisition
- Stress Analysis

NOTE

If WallMotion Segment Score is not displayed on the screen, select the “WallMotion” preset in the **Utility > Measure > M&A > Plot > Available Folders and Measurements**.

Getting started with a stress study

1. After selecting the appropriate application and probe and connecting ECG, press the **Protocol** tab on the Touch Panel. The protocol screen displays the layout of the default stress protocol for the current probe. This layout is also known as a template.


Table 6-2 Protocol Tab

Parameter	Description
Analyze	Display the Analysis screen
Template Editor	Display the template editor screen
Add Level	Add Level to the template
Delete Images	Delete the selected image
Move Image	Move the selected image to the another cell
Sync. Select	Synchronize the selected images.
End CC	End Continuous Capture
Begin/Cont.	Begin or continue the acquisition

Parameter	Description
Template	Display the template list
T1	Display/Hide the timer T1
T2	Display/Hide the timer T2
Cancel	Cancel Stress Echo


- To use the current template, press **Begin/Cont.** to initiate scanning.
To use another template, press **Template**. The template list displays.

Figure 6-4 Template List



- Trackball** to the desired template and press **Set**.
- The selected template displays.

Figure 6-5 Template (Example)



- a Level
- b Projection
- c Current Acquisition (green)

- Press **Begin/Cont.** to initiate scanning using the new template.


Image acquisition

Images are acquired in a pre-defined order, according to the selected template. The highlighted (green) cell of the template matrix, displayed in the Clipboard window, indicates which view is currently being acquired.

The names of the view and levels for the current cell are displayed on the top left corner of the image area and under the template matrix.

Acquisition Screen

Figure 6-6 Acquisition Screen



- 1 Current View Level
- 2 Timer
- 3 Template Matrix
- 4 Current View (Green cell)

Starting acquisition

1. Select the template.
2. Press **Begin/Cont.**
3. Perform a scan that conforms with the view that is highlighted in the template matrix on the Clipboard window.
4. Press the Store key.
 - If the actual stress level is configured to preview the cine loop before storing, use the cine loop controls to select the most appropriate heart cycle and if desired adjust the loop markers. Press **Store** again to save the selected cine loop.
 - or
 - If you do not want to store the cine loop, press **Freeze** to cancel. Return to the scan screen.

- If the actual stress level is not configured to preview the cine loop before storing, the system automatically stores the last heart cycle.

Stress levels can be configured for side-by-side display/comparison of the reference loop from baseline or previous level and the loop to acquire.

5. After storing the cine loop, the system automatically highlights the next view in the matrix to be acquired.
6. Repeat previous steps until all required views are completed.
7. If you select Auto Start Analysis on the Template Editor for this template, a dialogue asking “Do you want to start protocol analysis now?” displays when the last acquisition is complete. If you select Yes, the Stress Echo Analysis screen is displayed.

The template used can be configured so that analysis automatically starts by displaying the first protocol group. The wall segment scoring diagrams for each view is displayed in the Parameter window on the left side of the screen.

If the **Protocol** tab is selected during acquisition, the following Touch Panel displays.

Table 6-3 Protocol Tab during acquisition

Parameter	Description
Stop	Stop Stress Echo.
Pause	Pause Stress Echo. The template matrix continues in display. Even if you press P1, the cine loop does not store to the matrix.
Select Cycles	The Continuous Capture Selection screen is displayed (only available in Continuous Capture mode).
Analyze	Enter Analysis screen.
Template	Enter Template screen.
Add Level	Add level to the template.
T2	Display (Start)/Hide Timer T2.

Selecting a view during acquisition

A fixed protocol is provided for scanning, based on the selected template. The system automatically highlights the next view to be acquired in the template matrix, as images are stored. However, the order of scanning may be changed manually as follows:

Manual selection of a view during acquisition

1. Use the **Trackball** or the arrow keys on the alphanumeric keyboard to move the cursor to the cell that represents the view to be acquired.

The selected cell in the template matrix, highlighted in red, indicates the non-default position. When blinking, it contains a previously-stored acquisition.
2. Scan and save the selected loop as explained in the previous section.

After storage, the system automatically highlights the next available view to be acquired.

Moving an acquired image

An image can be moved from one cell to another during acquisition.

Procedure 1

1. When in the Protocol screen, press **Move Image**.
2. Use the **Trackball** to move the cursor to the desired image.
3. Press **Set**.
4. Use the **Trackball** to move the cursor to the destination cell.
5. Press **Set**. The image is moved from the source cell to the destination cell.

Procedure 2

1. In the Protocol screen, use the **Trackball** to move the cursor to the cell containing the image to move (source cell).
2. Press and hold down **Set**.
3. With the **Set** key still depressed, move the **Trackball** to the desired cell.
4. Release the **Set** key. The image is moved from the source cell to the destination cell.

If the destination cell contains an image, the images from the source and destination cells is exchanged when moving an acquired image.

Timers

Two timers can be displayed in the Stress mode acquisition screen, beside the template matrix.

Timers

- T1 displays the elapsed time from the start of stress examination.
- T2 starts when entering live scanning on the second stress level.

Both T1 and T2 timers can be manually stopped and restarted during the acquisition.

The display of T1 and T2 is user-configurable.

NOTE

If you activate the Timer in Stress Echo, the T1 timer is displayed in the lower left-hand corner of the image area after exiting Stress Echo.

Continuous Capture mode

Continuous Capture mode enables the user to perform acquisition continuously for all views at any level depending on the selected template configuration. Continuous Capture consists of temporary saving images acquired in a storage buffer. To enable best possible use of the limited storage buffer capacity, a Pause/Capture mode is provided, as opposed to the normal Freeze/Scan mode. The Pause mode enables scanning and live display on the screen, without any capture, thereby leaving the buffer available.

To run Continuous Capture, the user has to select a template where this feature is activated.

The buffer bar

When entering a level with Continuous Capture enabled, a buffer bar displays in the window. The Buffer bar displays the following information:

- The scanning state
 - Pause (live scanning without storing)
 - Capture (live scanning with storing to buffer)
- The percentage of the buffer that is filled
- The buffer filling progression showed by a filling gage
- The capturing sessions, reflected by the red lines along the buffer bar

Figure 6-7 Buffer Bar



Controlling the capture process

When entering a stress level with Continuous Capture enabled, the system is automatically set in Pause mode.

1. Press **Store** to start image capture.
“Capture” is displayed in the buffer bar, the gage starts filling and the percentage of filled memory buffer increases.
2. Press **Store** again to stop capture.
“Pause” is displayed in the buffer bar.

When 90% of the memory buffer is filled up, the text display in the buffer bar turns red.

The system enters Freeze mode automatically once the buffer is full and the captured loops display in the Continuous Capture selection screen.

Activating Continuous Capture

1. Do all your pre-stress acquisitions in the Cardiac application.
2. Press the **Protocol** tab to enter the Stress Echo mode. The Protocol screen displays.
3. Press **Template**. The template list displays.
4. Select the template **Exercise 2x4** from the list.
5. Press **Begin/Cont.**
6. Acquire the resting loops in all four views.

NOTE

Use the **Store** key to store the images.

7. Once the fourth loop is acquired, the system enters into a waiting mode where Continuous Capture is in a pause state awaiting the patient to exercise.

8. When the patient is back on the bed, press **Store**. The Continuous Capture acquisition starts.
9. Acquire all your views.

The memory buffer gage increases. When memory exceeds 90%, the percent number turns red.

10. Press **Freeze** to finish.


11. Press **Select Cycle**.

The Continuous Capture selection screen displays.

If the buffer is full, the system automatically displays the Continuous Capture Selection screen.

Refer to the next section if additional image acquisitions are necessary after the buffer is full.

Figure 6-8 Continuous Capture Selection Screen



12. Assign the cine loops to the four views.

a - Trackball to the desired loop.

- b** -Press **Set**. A drop-down menu appears with the available choices.

Figure 6-9 Drop-down Menu



- c** -**Trackball** to the appropriate view.

- d** -Press **Set**.


NOTE

To access additional cycles, use the arrow keys on the lower left portion of the select cycle screen.

- e** -Continue these steps until all views are selected.

- f** - Select **Done** when complete. A dialogue window displays, asking whether the entire Continuous Capture acquisition should be saved.

Figure 6-10 Dialogue Window



13. Press **Delete** to discard the loop or press **Store all** to keep the entire loop.
14. Perform Analysis and scoring.

Continuous Capture with additional image acquisition

If the buffer is filled up before all the image acquisitions are done, additional loops can be stored in the clipboard before doing image assignment to the views.

1. Perform the Continuous Capture. *Activating Continuous Capture* on page 366 (Steps 1 to 11).
2. In the Continuous Capture selection screen, press **Select Later**.

The Continuous Capture screen displays.

3. Perform the additional acquisition.
4. In order to resume the stress echo exam and assign loops for the views from the Continuous Capture buffer, press **Protocol**. If not displayed, select the template **Exercise 2x4** from the template list.
5. Click the continuous capture images on the Protocol Template screen.

The Continuous Capture selection screen displays.

6. Assign the cine loops to the view. *Activating Continuous Capture* on page 366 (Step 12 a - f).
7. Press **Delete** to discard the loop or press **Store all** to keep the entire loop.

The normal procedure is to discard the loop. The loop is very big and requires a lot of disk space.

8. Perform Analysis and Scoring.

Postponed image assignment

The assignment of the cine loops to the view can be done on a later stage on a stored Continuous Capture acquisition.

1. Perform the Continuous Capture. *Activating Continuous Capture* on page 366 (Steps 1 to 11).
2. Press **Store all**.

The entire Continuous Capture acquisition is stored. The examination can be ended and the image assignment, analysis and scoring can be done later.

3. Re-open the examination, if necessary.
4. Press **Protocol**. The Protocol screen displays.
5. Click the continuous capture images on the Protocol Template screen.

The Continuous Capture selection screen displays.

6. Assign the cine loops to the view. *Activating Continuous Capture* on page 366 (Step 12 a - f).
7. Select **Done**.
8. Perform analysis and scoring.

Restart capture from the Continuous Capture Selection

Press **Restart Capture**.

The recording in memory is deleted and the Continuous Capture starts again.

Resume Continuous Capture


Press **Continue Capture**.

Resumes Continuous Capture recording (only if the Continuous Capture buffer is not full).

Assigning and storing the cine loop

The cine loops captured in the buffer are assigned to the stress protocol views and stored from the Continuous Capture selection screen.

Figure 6-11 Continuous Capture Selection Screen




- a. Highlighted loop
- b. Cycle number and total number of cycles
- c. Blue Gauge: Position of the highlighted loop within buffer area.
- d. Navigation Controls: << < > >> (back to first selection, back to previous selection, forward to next selection, and forward to final selection).

Assigning a cine loop to a view

1. Use the **Trackball** to move the cursor to the desired cine loop in order to assign it to a particular view of the stress template.
The frame of the loop is highlighted.
2. Press **Set**.

A pop-up menu displays with the view names of the template.

Figure 6-12 Loop Assignment



- a. Already assigned view
- b. Views pop-up menu

NOTE

A checkmark appears on the Views pop-up after you have assigned a view to an image.

3. Use the **Trackball** to select the required view name.
4. Press **Set**.
The name of the view displays above the timers in the window.
5. Repeat steps 1 through 4 to assign loops to the other views of the level.
6. Press **Done** when complete.

A dialogue window displays asking whether the entire Continuous Capture acquisition should be saved.

7. Press **Delete** to discard the loop or press **Store all** to keep the entire loop.

The normal procedure is to discard the loop. The loop is very big and requires a lot of disk space.

Post Acquisition Features

Post acquisition, you can utilize Raw Data to adjust the following in B-Mode:

- Zoom

Advanced Features

- SRI HD
- Rejection
- Frame Average
- TGC
- Maps
- Dynamic Range
- Gain
- Rotation

You can also take measurements post Stress Echo acquisition.

Analysis

Analysis consists of viewing previously saved loops and assigning scores to each cardiac segment, in order to quantify the function of the muscle or wall segment.

Depending on the protocol configuration, the analysis stage can start manually or automatically after completion of the stress test. In this case, the usual procedure consists of sequentially opening all image groups (if defined) and performing scoring from image to image.

The quad screen is the standard display for comparing heart cycles. The heart cycle loops in the display are synchronized to enable comparison. Each loop in the quad screen can be magnified, using the zoom control.

Image Selection for Analysis

Images can be selected manually or from a pre-defined group in the Protocol screen.

Selection of Images from a group

If groups of images have been defined in the protocol template, you can select a group of images for analysis and sequentially analyze all images from all groups from within the Analysis screen.

1. In a stress examination, press **Protocol**. A preview of the acquisition displays.
2. Press **Analysis**. A pre-defined group appears in the display with a Wall Segment window on the left.
3. To advance to other groups, use the **Trackball** to move the cursor to the arrows at the bottom of the Wall Segment window. Select an arrow to advance to another group. For further clarification, see callout F in *Figure 6-13 Analysis Screen* on page 373.

Manual selection of images from Analysis screen

1. When currently in the protocol analysis screen in the Stress analysis quad screen, hold down the **Shift** key while performing Steps 2 through 4.
2. Use the **Trackball** to move the cursor to the first image to select in the template matrix.
3. Press **Set**. The frame of the selected loop is in the Stress analysis screen and the next window in the quad screen is automatically selected.
4. Repeat step 2 and 3 to select other images.

- Depress **Shift**.

Manual selection of images in the Protocol screen

- In a stress examination, press **Protocol**. A preview of the acquisition displays.
- Use the **Trackball** to move the cursor to the first image to select.
- Press **Set**. The frame of the selected loop highlights.
- Repeat Steps 2 and 3 to select other images.
- Press **Analyze** to open images in the Analysis screen.

Scoring acquired loops

- After image selection, press **Analyze**.
The Stress Echo analyze screen displays.

Figure 6-13 Analysis Screen



- Wall segment diagram
 - Selected loop (Highlighted frame)
 - Displayed loops (Highlighted frames)
 - Exit Wall motion scoring
 - Change page or enter next image group
- Use the **Trackball** to move the cursor to a segment in one of the scoring diagrams and press **Set**.
The Score pop-up list displays.
 - Use the **Trackball** to move the cursor to a score.
 - Press **Set**.
The score displays in the relevant segment area in the diagram.

NOTE

To edit a score, select it and choose a new score.

5. Repeat step 1 through 3 to score relevant segments.
6. Press the **Change Page** arrow to display the next group of images.
7. Repeat step 1 through 3 to score relevant segments on the new loops.

Figure 6-14 Segment Scoring



NOTE

Since Cine changes into sync mode, subsequent scans are also synchronized. Exit sync mode from the Cine menu.

Editing/Creating template

The stress package provides protocol templates for exercise as well as pharmacological stress examinations.

The user can create new templates or modify existing templates to suit the individual needs. Up to ten projections and fourteen stress levels can be created in a template.

Templates created may be temporary, used only during the current examination, or saved as new templates, for future use and reference.

The editions that may be performed include:

- Adding/Deleting levels and projections.
- Assigning new labels to levels and projections.
- Defining level options.
- Defining new groups.

Templates are edited/created from the Template editor screen.

Entering the Template editor screen

1. Press **Protocol** to enter the stress echo mode.
2. Press **Template**. The template list displays.
3. Use the **Trackball** to select the Template Editor.
4. Press **Set**. The Template Editor screen displays.

OR

1. Press **Protocol** to enter the stress echo mode.
2. Press **Template Editor** on the Touch Panel. The Template Editor screen displays.

Template Editor screen overview

Figure 6-15 Template Editor Screen




Table 6-4 Template

Parameter	Description
Template	Select a pre-defined template from the pull-down menu. The protocol template preview updates accordingly.

Table 6-5 Protocol Template Preview

Parameter	Description
Protocol Template Preview	- Displays an updated preview of the template accordingly to the settings applied. - To change Projection and Stress level labels, select a pre-defined label from the pull-down menu or press Set in the actual label field and type a new name.

Table 6-6 Template Settings

Parameter	Description
Template Settings	- Cycles: select the number of cine loop heart cycles to store for each level from the pull-down menu or enter the desired value manually. - Continuous Capture: Checking this parameter enables continuous image acquisition throughout the level. The images acquired are temporarily stored in the system's storage buffer. - Preview of store: Checking this parameter enables review and adjustment of cine loops before store. - Show reference: Checking this parameter displays a dual screen with the reference level (first or previous level) on the left and the live image on the right.

Table 6-7 Scan Mode

Parameter	Description
Scan Modes	2D, Color, PW (Pulsed Wave Doppler), CW (Continuous Wave Doppler), MM (M-Mode), Color MM, Color PW, Color CW

Table 6-8 Other options

Parameter	Description
Other Options	- Grid Size: Enter the number of levels and projections for the selected template. - Timers: If you check this parameter, starts T1 and T2 timers automatically. - Auto-start analysis: If you check this parameter, displays the Stress Echo Analysis when the last acquisition is performed. - Show Systole in Analysis: When selected, the systolic part of the cardiac or ECG cycle is only displayed. The whole cycle is not displayed. - Smart Stress: Check Smart Stress to store a subset of the image acquisition settings (e.g., Zoom, Gain, Compress, etc.) for each view in the protocol. Smart Stress enables you to set image acquisition settings for each view at a baseline level and automatically get the same image settings in the corresponding views in the next levels. - Reference image: When Show Reference is selected, selects either corresponding baseline loop or corresponding loop from the previous level to be displayed as reference image during acquisition.

Table 6-9 Pre-defined groups

Parameter	Description
Pre-defined groups	- Shows the image groups created. - New group: Creates a new image group. Select the desired images on the template preview. - Update group: Edits a selected group after new loop selection on the template preview. - Delete group: Deletes a selected group.

Selecting a base template to edit

1. Select the base template from the template pull-down menu on the upper left corner.
2. Press **Set**.
The selected template displays in the protocol template preview field, showing the levels, projections and their labels.

Adding/Deleting levels and projections

1. Enter the number of levels and projections in the Grid size field.
The new grid size displays in the protocol template preview field.
2. Press **New Template** to create a new template.
or
Press **Save Template** to update the base template.

Display timers

Check the box(es) to display timer(s) as specified.

NOTE

The timers can also be started or stopped at any time during stress examination by using the T1 and T2 Touch panel key.

Start analysis automatically

Check **Auto Start Analysis** to display the Stress Echo Analysis screen when the last acquisition is performed.

Smart Stress

Check Smart Stress to store and automatically reuse a subset of the image acquisition settings in the baseline level view in the corresponding views in the next levels.

Configuring levels

The following options can be set up for each level:

Number of cycles to be stored in the cine loop:

Enter or select the desired number in the Cycles field.

Continuous Capture

Check Continuous Capture if continuous image acquisition throughout the level is desired.

When Continuous Capture is selected, preview of the cine loop and reference display during acquisition are not possible.

Preview of store

Check Preview of store if review and adjustment of cine loops before storage is desired.

Show reference

Check Show reference if the display of the corresponding reference loop is desired during acquisition (Dual screen mode).

Adding a group

1. In the Protocol template preview field, select the cells to be part of the group.
2. In the Pre-defined group field, press **New group**.
A dialogue box displays to ask the user to enter a name for the new group.
3. Enter the group name.
4. Press **OK**. The new group displays in the pre-defined group field.

Updating an existing group

1. In the Pre-defined group field, select the group to edit.
The selected cells are highlighted in the Protocol template preview field.

NOTE

The selected group is highlighted by a Green frame.

2. Either select a new cell(s) to add to the group or deselect an existing cell(s) to remove from the group.
3. Press **Update group**.
The display in the Protocol template preview field is updated accordingly.

Deleting a group

1. In the Pre-defined group field, select the group to delete.

NOTE

The selected group is highlighted by a Green frame.

2. Press **Delete group**.
The group is removed from the list in the pre-defined group field.

Specifying Scan Mode for each Projection

Specify the Scan Mode for each Projection: 2D (B-Mode), Color Flow Mode, M-Mode, Color M-Mode, PW Mode, Color PW Mode, CW Mode, or Color CW Mode.

Saving the Template

You can save the template using controls at the bottom of the Template Editor page, or use the controls on the Touch Panel.

Table 6-10 Template Editor Saving Options

Parameter	Value
New Template	Select this option to create an entirely new template.
Save As Template	If you would like to create a new template based on the existing template with your modifications, select to Save this Template As, and give it a name.
Save Template	Select this option to save the default template with your modifications.
Delete Template	Select this option to delete a template.

Wall Motion Segment Setup

You can set up the following parameters for Wall Motion Segment in the Utility screen (**Utility** > **Measure** > **Advanced** > **Cardiac**).

Figure 6-16 Cardiac Wall Motion Segment Setup




Table 6-11 Wall Motion Segment Parameters

Parameter	Value
WMS freeze loop at ES	Specify to freeze the Loop at End Systole
WMS Segment Model	Select 16 or 18 segments
WMS initial scoring	Undefined or Normal
WMS scoring legend	ASE, ASIA or European

Utility Application Settings for Protocol

Figure 6-17 Protocol Setup

The screenshot shows the 'Settings' tab of the 'Utility Application Settings for Protocol' window. The 'Preset' is set to 'Abd'. Under 'Image Control & Display', 'Join Dual Image for Linear' is checked. Under 'When Entering Dual Image...', both 'Duplicate Frozen Image to Opposite Side' and 'Duplicate Live Image to Opposite Side' are unchecked. Under 'Patient Info', 'Title Bar Line 1' is 'Last,First Name', 'Title Bar Line 2' is 'ID', and 'Title Bar Line 3' is 'Empty'. Under 'Comments', 'Active Function at Freeze' is 'None'. Under 'Footswitch', 'Left' is 'Print 1', 'Middle' is 'Freeze', and 'Right' is 'Print 3'. The 'Protocol' section is highlighted with a blue box; it contains 'Show Protocol Tab' (unchecked) and 'Template' set to 'Pharmacological 4x4'. Under 'ECG', 'Show ECG Tab' is unchecked. Under 'ELASTO', 'Show Quality Bar' is checked and 'Show Quality Graph' is unchecked.

Table 6-12 Protocol Parameters

Parameter	Description
Show Protocol Tab	Show/Hide the Protocol tab for that preset (Bicycle Normal, Bicycle Sporty, Contrast Pharmacological, Pharmacological 4x4, Pharmacological 8x5, Exercise 2x4, Exercise 2x4 B, Pharmacological US 4x4, or User-configured).
Template	Select the default template.


Reports

If you set up the Wall Motion Analysis field on the Report, you can insert the results.

Advanced Features

Select Report to view either the Bull's Eye and/or Cut Plane Report.

Figure 6-18 Bull's Eye and Cut Planes Report Sample



ECG Option



WARNING

DO NOT USE the physiological traces of the Ultrasound system for diagnosis and monitoring.



CAUTION

Use GE HealthCare provided ECG cable to ensure the protection performance against the effects of a discharge of a cardiac defibrillator to the patient.



CAUTION

The operating environment requirement of the ECG is inconsistent with the console, and the operating environment temperature with the ECG should be controlled above 10°C.



CAUTION

- The conductive parts of the electrodes used in the applied part and their connectors (including the neutral electrode) should not come into contact with any other conductive part, including ground.
- After the defibrillator stimulates the patient, the ECG needs 5 seconds for recovery.
- The quality of the ECG trace depends on the stability and conductivity of the electrodes during the test, especially during high stages when the patient's movements can cause artifacts.
- Make sure that the lead wires do not swing.
- The device is not waterproof. Do not expose the device to water or any kind of liquid. Maintain in a dry place:
- Worn or damaged patient cables are the most common cause of poor ECG signals. ECG signals (or wave patterns) that consistently contain noise or artifact may suggest need for ECG wire or cable replacement.
- Store the device in a dry place.
- Always protect the recorder from coming into contact with moisture. In rain or snow conditions, protect the recorder from bad weather elements by wearing the recorder inside a coat.

Intended Use for ECG

During the echocardiography examination, the electrocardiogram (ECG) is connected to record the electrical activity of the heart synchronously, so as to analyze the regular changes of the phase of cardiac activity.

ECG connecting

The ECG module is connected to the system via USB ports.

Figure 6-19 Connecting ECG to the system



CAUTION

Install ECG Hardware Option before connecting the ECG trunk cable to the system.

NOTE

The conductive parts of the electrodes used in the applied part and their connectors (including neutral electrodes) should not touch any other conductive parts, including the earth.

Restart the system after ECG module is connected.

To activate ECG Option


Go to **Utility > System > User Configurable Key** to define a key for activating ECG.

To display the ECG Signal on the monitor, go to **Utility > Imaging > General > ECG > ECG Display**.

Operation


1. Connect the standard A-type plug of the USB cable (Detail 1) to the USB connector of the host.

Figure 6-20 USB cable




2. Connect the Mini-B type USB plug of the USB cable (Detail 2) to the input of the ECG-USB1 acquisition device (Detail 1).

Figure 6-21 ECG-USB1



3. Connect plug of detachable patient cable (Detail 3) to the round connector on ECG-USB1 (Detail 2).

Figure 6-22 Detachable patient cable



4. Apply the leads to the patient (Right Arm, Left Arm and Right Leg) according to their markings (R(RA), L(LA) and N(RL) respectively).
5. Operate through the host system.

Patient preparation

The ECG traces quality depends very much on the stability and conductivity of the electrodes during the test, especially during high stages when the patient movements can cause artifacts. Here are some basic rules to ensure good electrical contact:

- Shave hair at the electrode contact points
- Use high quality liquid gel electrodes
- Make sure that the lead wires do not swing

Attach the leads as shown in *Figure 6-22 Detachable patient cable* on page 386. Detachable 3 leads each 100cm with clip connectors on the patient end designed for use with disposable ECG electrodes.

Cleaning and Disinfection



CAUTION

ECG is not waterproof. Do not expose the device to water or any kind of liquid. Maintain in a dry place.

**CAUTION**

Before performing any cleaning and disinfection, always ensure that the product is disconnected from the patient and other equipment.

To avoid the risk of damaging the product, do not use any automated cleaning and disinfection procedures.

Only use approved cleaning and disinfection agents.

Frequency of visual inspections, cleaning, and disinfection

- ECG modules should be cleaned weekly - ECG lead wires should be cleaned and disinfected after each patient use.
- Visually inspect products before and after cleaning and disinfection procedures.
- Damaged products should not be used. Evidence of product damage may include, but is not limited to, corrosion, discoloration, excessive scratching, flaking, cracking, and wear.

Permitted Cleaners and Disinfectants

PDI Super Sani-Cloth Disposable Sterilization Swab

Cleaning

1. Put on new gloves.
2. Perform a pre-cleaning and wipe the surface of the product to ensure that all surfaces are cleaned evenly. Repeat cleaning until the product is visibly clean.
3. Use cotton swabs or other suitable cleaning accessories (such as soft brushes) to clean hard-to-reach areas. Do not use any sharp tools to clean the product.
4. The product can be rinsed manually thoroughly with a clean, soft, lint-free cloth saturated with warm tap water (water temperature range 0 to 50°C, 32 to 122°F).
5. Allow the product to air dry until visibly dry or wipe it thoroughly dry with a clean, soft, lint-free cloth. Drying time may vary depending on environmental conditions.

Disinfection

1. Put on new gloves.
2. Wipe the entire surface of the product. The treated surface must remain visibly wet for at least two minutes (Super Sani-Cloth) or according to the disinfectant manufacturer's instructions.
3. Use cotton swabs or other suitable cleaning attachments (e.g. soft brushes) to clean hard-to-reach areas. Do not use any sharp tools to clean the product.
4. The product can be thoroughly rinsed manually using a clean, soft, lint-free cloth saturated with warm tap water (water temperature range 0 to 50°C, 32 to 122°F).
5. Allow the product to air dry until visibly dry or wipe it thoroughly dry with a clean, soft, lint-free cloth. Drying time may vary depending on environmental conditions.

Troubleshooting - Noisy ECG Signal

"Noise" refers to any degradation of the ECG signal that makes it difficult to accurately detect and classify beats. Causes of noise, such as artifact and electrical interference, should be avoided whenever possible.

Some common causes of noisy ECG signals include:

- Poor skin prep.
- Dried electrode gel.
- Muscle artifact caused by shivering, movement or tremors.
- Baseline wander caused by excessive chest movement, or the electrical differences between two brands of electrodes.
- Respiration artifact caused by thoracic or abdominal movement of both spontaneous and ventilated breathing patterns.
- Nearby electrical equipment.

Table 6-13 Noisy ECG Signal


Problem	Appearance	Cause	Corrective Action
Baseline Wander	Rhythmic Up-and-down movement of ECG baseline	Movement of the Patient. Improperly applied electrodes.	Make sure patient is comfortable and still Reapply electrodes
Irregular Baseline	Rough, jagged baseline	Poor electrical contact. Faulty or dry electrodes.	Reapply electrodes, using proper technique. Check for loose connections on lead/cable Apply new electrodes
Muscle Artifact	Fuzzy, irregular baseline	Tense, uncomfortable patient. Tremors, Diaphoresis. Poor electrode placement.	Make sure patient is comfortable and still. Check that electrodes are applied on flat, on-muscular areas of the torso. Reapply electrodes if necessary.
Poor Electrical Contact	Trace switching from high to low in steps and/or dashed trace	Loose electrodes Defective leads/cables	Change all electrodes, using good skin prep. Replace leads/cables.
Power Line Interference (50/ 60 Hz)	Regular saw tooth baseline.	Poor electrode placement. Possible non-grounded instrument near patient.	Reapply electrodes. Check grounding of equipment near patient.

ECG Trace Monitor Display

The scanned image is synchronized with the ECG trace. In Doppler or M-Mode, the traces are synchronized with that particular mode's sweep.

The user can control the gain, position and sweep speed of the traces.

Figure 6-23 ECG Trace Display



1. ECG
2. Auto Heart Rate Display

ECG menu

The ECG menu provides for control of the input signals.

Without the ECG option, the ECG menu will not be displayed.

NOTE

After connect ECG module to the ultrasound system, please restart to activate the ECG option.

Figure 6-24 ECG Touch Panel




Table 6-14 ECG parameters

Parameter	Description
Sweep Speed	Change the speed of the trace. The sweep speed of the physio signal on the B-Mode image can be set independent of the timeline (Doppler and M-Mode) sweep speed.
ECG Trigger	Enables intermittent imaging based on the ECG. The trigger location(s) relative to the R trigger are set with the Delay Time key. Press ECG Trigger and select one of the options (None, Trig1, Trig2, and Both) and adjust the delay time using the Delay Time key. ECG Trig 1 specifies the delay (ms) from R-wave to triggered frame. ECG Trig 2 specifies the delay from R-wave to second frame. Both activates ECG Trig 1 and ECG Trig 2 simultaneously. Trig 2 must be greater than Trig 1 for dual triggering (Both) to be active.
ECG Display	Provides the ability to turn on the ECG trace and Auto Heart Rate for display on the monitor.
ECG Trigger Period	The control specifies the number of heart cycles (R-waves) that are skipped between ECG triggers. The default is 1 or no skipping.

Parameter	Description
Delay Time	In ECG Trigger Mode: If only ECG Trig1 or ECG Trig2 is selected via the ECG Trigger key, the Delay Time key controls the R-Delay time of the active trigger. If both triggers are selected (Both), press this key to toggle ECG Trig1 and ECG Trig2 and rotate the key to change the delay time. Once the trigger is set, the snap shot image is displayed each time the update line passes the active trigger(s). In Timer Trigger Mode: Rotating the knob changes the delay time between images.
ECG Gain/ Position	Allows for the amplitude control of the ECG trace or allows for the vertical positioning of the ECG trace on the image display. Press the knob to toggle between Gain and Position. The default is Gain.

Using 3D (Option)



WARNING

DO NOT scan any pacemaker patient using the sensor device. The magnetic fields emitted from the device may interfere with the pacemaker operation.

Easy 3D is compatible with every 2D transducer using a freehand acquisition to generate a volume dataset.

Table 6-15 3D Package Options

3D Type	Description	Sensor/No Sensor	Available Tabs
Easy 3D	Designed for rendering B Mode images, e.g., Baby Face scans.	No sensor	3D Acquisition, Easy 3D, Movie
Advanced 3D	Designed for rendering B Mode images, e.g., vessel trees.	No sensor	3D Acquisition, Easy 3D, Advanced 3D, Movie

3D Volume datasets are allowing the navigation in the 3D cube itself and providing access to the 3 different main planes - axial, sagittal and coronal.

Easy 3D

Acquiring a 3D Scan

To acquire a 3D scan,

1. Optimize the B-Mode image. Ensure even gel coverage.
2. Press the **4D/3D** key on the control panel.

NOTE

Set appropriate values for PreAcq and Scan Plane. Also, set the scan distance before scanning.

3. To start acquiring the image, press **Start** key.
4. To perform a parallel scan, scan evenly. To perform a sweep (fan) scan, rock the probe once. Note the distance of the scan.
5. The 3D volume of interest (VOI) is dynamically assembled on the right side of the screen.

NOTE

If the image stops before you're done scanning, start acquiring the 3D volume of interest again.

6. To complete the 3D scan, press **PreMode** key.

NOTE

You can also press Freeze, but then you need to also press the 3D key to obtain the final render.

3D Notes

- Adjust the 3D data set brightness with B-Mode Gain.
- Use Colorize to change the color of the active data set.
- Use Zoom to increase the zoom factor of the active data set.
- Vertical lines may be seen in a resliced image. This usually happens when you scan too fast or if the scan distance is set to a high value.

Scan more slowly, adjust the frame rate for a faster rate or adjust the scan distance.

Manipulating the Volume of Interest


Imagine you are able to manipulate the 3D volume of interest (VOI) in your hand. The 3D VOI is a tangible anatomical object that you can see and manipulate easily using the Trackball and Set control panel keys.

Practice positioning the pointer at different places within the 3D VOI. Highlight different colors (white, red, yellow, or green). Press Set to select a VOI for manipulation. Use the hand to manipulate the 3D VOI.

Rotating the 3D VOI Left/Right or Forward/Backward

You can rotate it left to right or right to left. You can rotate it forward/backward. Press right **Set** key when the white pointer finger is positioned on the white box. Move the closed white hand to manipulate the 3D VOI.


Figure 6-25 Manipulating the 3D Volume of Interest (White Hand)



Moving Through the 3D VOI

You can move through the 3D VOI using the red hand. Press Set when the red pointer finger is positioned on the red box. Move the closed red hand to move through the 3D VOI.

Figure 6-26 Moving through a 3D Volume of Interest (Red Hand)



NOTE

Any plane in the volume can be made active (highlighted with red box) by clicking on it.


Viewing Specific Portions of the Anatomy

You can pull back tissue to view specific portions of anatomy using the yellow hand. Press Set when the yellow pointer finger is positioned on the yellow box. Move the closed yellow hand to manipulate the 3D VOI.

NOTE

This actually moves an edge. A yellow hand appears only when the pointer is on an edge of the VOI.

Figure 6-27 Manipulating the Edge of a 3D Volume of Interest (Yellow Hand)


**Pulling Back a Corner of the VOI to View Specific Anatomy**

You can pull back a corner to view specific portions of anatomy using the green hand. Press Set when the green pointer finger is positioned on the green box. Move the closed green hand to manipulate the 3D VOI.

NOTE

A green hand appears only when the pointer is on a corner of the VOI.

Figure 6-28 Manipulating a Corner of the 3D Volume of Interest (Green Hand)



3D Acquisition Parameter Description

Table 6-16 3D Acquisition Description and Instructions for Use

3D Parameter	Description
Application Presets	Selections: None, OB - Baby Face None. No application preset applied. OB - Baby Face. After having scanned in this mode, certain rendering parameters are set automatically. The gray surface mode is activated and the texture mode is switched off. The gray surface mode values for opacity and threshold are set automatically according to the datasets histogram. User1, User2, User3. User defined preset.

3D Parameter	Description
Acquisition Mode	Selections: Sensorless Parallel, Sensorless Sweep Sensorless Parallel. In this mode the probe must be moved during 3D data acquisition without angling it. You should scan the object you want to render in 2-4 seconds. The speed at which you scan should be constant. No sensor is mounted on the probe. Since the time for post-processing depends on the acquired number of frames, it is recommended that you check the frame rate. Low frame rates result in fewer acquired frames for the 3D dataset which results in intensive post-processing (interpolation). Therefore, low frame rate = long post-processing. Sensorless Sweep. In this mode the probe must be moved to a position where you can clearly see a middle cut of the object you want to scan and render. Tilt the probe to about 30 degrees until the object you want to scan disappears. Start the acquisition and tilt the probe over a distance of around 60 degrees until the object disappears again. The entire scan time should be around 2-4 seconds. During the sweep, the probe may not be moved parallel, just tilted. No sensor is mounted on the probe. Before starting an acquisition, take care that the transmitter is positioned correctly during data acquisition and that the transmitter cannot move.
Scan Plane	Selections: Front to Back, Side to Side Front to Back. After having scanned in this mode, the rendered dataset is shown in a frontal view. For acquiring a fetal face in sagittal cuts, use this mode. Side to Side. After having scanned in this mode, the rendered dataset is shown from a side view. For acquiring a fetal face in coronal cuts, use this mode.
3D	Starts the rendering process.
Scan Distance	Adjusts the distance covered during the scan. Depending on the real width of a scan acquired during a sensorless 3D acquisition, the volume of interest's width can be enlarged or reduced. You can adapt the form of a fetal face if the baby's head looks oval instead of round. The assumed default width of a parallel scan is 6 cm; of a sweep scan 60 degrees.
Display Format	Change the layout between 50/50 and only 2D.

NOTE

Selection of user preset shall be effective only while 3D mode is active. Exiting 3D mode and activating again then 3D preset shall be changed to the default, even without New Patient or changing application.


NOTE

When 3D image is recalled, no 3D preset shall be active and parameters are recalled from image file.

NOTE

Default Scan Distance, Opacity and Threshold may not consistent and change per scan. After User Preset was saved and recalled, Opacity and Threshold are consistent.

Figure 6-29 Easy 3D



Descriptions and instructions for using Easy 3D follow:

Table 6-17 Easy 3D Description and Instructions for Use

3D Parameter	Description
Reset	Resets the 3D volume of interest back to its original orientation.
Utilities	Select Average Off, Average Light, Average Medium, or Average Strong.
Undo	Undoes any manipulation you have done to your 3D dataset.
Auto Movie	Initializes the calculation and display of a 3D movie. A rotation of 30 degrees left and right around the actual image position (either the default position after acquisition or the position that was manually defined by manipulating the 3D volume of interest) is shown. For this 60 degree rotation, eleven images in steps of 6 degrees are calculated.
Scalpel	Structures, for example a part of the placenta hiding the view to a fetal face, can be cut out in a rendered image. All visible structures can be cut out. The option of 'erase inside' deletes all structures inside the marked region. The option of 'erase outside' deletes all structures outside the marked region. The region in the rendered image is marked with the right Set key. To define the contour of the region, press the right Set key for each vertex. To close the contour, double click the right Set key. As long as a contour is not closed, it can be traced back with the left Set key. The cut out process can be undone by the Undo Last function. As soon as the Apply button is pressed, a new dataset is generated.
Gray Surface	Activates the gray surface rendering mode. It leads to a transparent appearance of the object, generated by displaying only a surrounding shell of structures.

3D Parameter	Description
Texture	Activates the texture or photorealistic rendering mode. It creates a photorealistic appearance of the object. The shading depends on the orientation of the surface of the object. If both Texture and Gray Surface mode are switched on, the mixture percentage of both modes can be defined.
Render	Changes between the rendered image view and the view of a volume of interest. The volume of interest shows the acquired ultrasound images transformed into an isotropic rectangular coordinate system. The volume of interest can be manipulated as described above.
Threshold/Opacity	Threshold defines which gray values are used for rendering and which are considered noise. Opacity defines how strict Threshold is used for discrimination. A low opacity value creates a firmer appearance of the surface. A high opacity value leads to a transparent appearance of the rendered image.
Colorize/Contrast	Colorizes the 3D render or adds contrast to the 3D rendered image.
Overwrite User1/ User2/User3	Overwrite the application preset file with the changes you just made.
Orientation Marker	You can now specify/define, then add the following orientation markers while in 3D via the Orientation Marker key: • TRV Sup to Inf Ant Scan Prb Rt • TRV Inf to Sup Ant Scan Prb Rt • SAG Lt to Rt Ant Scan Prb Sup • SAG Rt to Lft Ant Scan Prb Sup • Defined Superior Inferior Left Right Anterior Posterior Cancel • None

3D Parameter	Description
Active Data	Manipulations of rendering parameters only have an effect on the data defined as Active Data. After having selected Active Data, a list of data is displayed, Gray Data or Inversion. Choose the data to be manipulated. Active Data is only available when you select both Inversion and Gray Data in Visible Data. Inversion Mode is only available for Black-and-White mode.
Visible Data	After selecting Visible Data, a list of data is displayed, Gray Data or Inversion. Choose the data you want to display. For example, if only Inversion is chosen, the B-Mode image is switched off in the rendered image and only inversion mode is displayed.

Advanced 3D

Descriptions and instructions for using Advanced 3D follow:

Table 6-18 Advanced 3D Descriptions and Instructions for Use

3D Parameter	Description
Reset	Resets the 3D volume of interest back to its original orientation.
Contrast	Adjust to increase or decrease the contrast of the image.
Colorize	Colorizes the gray scale image to enhance the eye's discrimination capability.
Utilities	Use smoothed volume for rendering the 3D. Strong = Most Smoothing.
Auto Movie	Initializes the calculation and display of a 3D movie. A rotation of 30 degrees left and right around the actual image position (either the default position after acquisition or the position that was manually defined by manipulating the 3D volume of interest) is shown. For this 60 degree rotation, eleven images in steps of 6 degrees are calculated.
Save Preset	Save as a user preset (User1, 2, or 3).
Undo	Undoes any manipulation you have done to your 3D dataset.
Tile	The display can be divided into 1, 2, 4, or 6 windows. Switching to a lower number of windows keeps the images from left to right.
Scalpel	Structures, for example a part of the placenta hiding the view to a fetal face, can be cut out in a rendered image. All visible structures can be cut out. The option of 'erase in' deletes all structures inside the marked region. The option of 'erase out' deletes all structures outside the marked region. The option of 'erase outside' deletes all structures outside the marked region. The region in the rendered image is marked with the Enter key. To define the contour of the region, press the right Set key for each vertex. To close the contour, double click the Enter key. As long as a contour is not closed, it can be traced back with the left Set key. The cut out process can be undone by the Undo Last function. As soon as the Apply button is pressed, a new dataset is generated and the cut is permanent

3D Parameter	Description
Active Data	Manipulations of rendering parameters only have an effect on the data defined as Active Data. After having selected Active Data, a list of data is displayed, Gray Data or Inversion. Choose the data to be manipulated. Active Data is only available when you select multiple items in Visible Data. NOTE Inversion Mode is only available for Black-and-White mode.
Visible Data	After selecting Visible Data, a list of data is displayed, Gray Data or Inversion. Choose the data you want to display. For example, if only Inversion is chosen, the B-Mode image is switched off in the rendered image and only inversion mode is displayed.
Define Axis	For certain display and measurement modes (Angular Plane Mode, Angular Volume Measurement Mode), an axis in the volume of interest is required. To define the axis, set the start point by using the Trackball to position one end of the axis and pressing the Enter key, then positioning the other end of the axis and pressing the Enter key.
Group Planes	Selections: Off, Main, Parallel Off. A VOI or a rendered image is displayed. The Render button changes between the rendered image view and the view of the VOI. The VOI shows the acquired Ultrasound images transformed into an isotropic rectangular coordinate system. Main. Three orthogonal cuts (with colored frames) of the acquired VOI are displayed after pressing Main. The VOI shows the acquired ultrasound images transformed into an isotropic rectangular coordinate system. On the left top of the image a complete VOI is displayed. It shows the position of the three orthogonal planes in the VOI. A green point displayed in each plane defines the point of intersection of the three planes. This point can be set to different positions in the planes by double clicking the Enter key. A plane can be moved parallel in the VOI by pressing the Enter key on the position of the green point and moving the Trackball up and down inside the plane. Parallel. In this mode all displayed VOIs get the orientation of the last modified volume. Normally four VOIs are displayed. It is possible to display six VOIs by increasing the number of displayed volumes in the Tile area. Between the first and the last VOI, the selected planes are parallel and equidistant. A modification on the plane in one VOI results in a parallel modification of the planes in all other VOIs. To move these planes, press and hold down the Enter key while moving the Trackball.

3D Parameter	Description
Type 1/2	Defines the rendering modes. Selections: Gray Surface, Texture, Maximum Intensity, Minimum Intensity, Average Intensity, and None. If both Type 1 and Type 2 rendering modes are switched on, the mixture of both modes can be defined. Gray Surface. Activates the gray surface rendering mode. It leads to a opaque appearance of the object, generated by displaying only a surrounding shell of anatomical structures. Adjust Threshold and Opacity as well. Texture. Activates the texture of photorealistic rendering mode. It creates a photorealistic appearance of the object. The shading depends on the orientation of the surface of the object. Adjust Threshold and Opacity as well. Maximum Intensity. Transparent appearance of the object. Generated by displaying the maximum gray values in the VOI. Minimum Intensity. The rendered image is generated by displaying the lowest gray values in the VOI that exceed the defined threshold. Dark anatomical structures, like cysts, can be shown in this mode. Average Intensity. Transparent appearance of the object. Generated by a summation of the gray values. None for Type 2. No second rendering mode is used in addition to the Type 1 rendering mode.
3DLandscape	Shows a combination of 2D slices and a 3D rendered image. After a color acquisition you can combine the 2D B-Mode image slices with a 3D rendered color image. This mode allows stepping through the B-Mode images along a vessel structure. The 2D slice can be moved with the Enter key. The Trackball symbol has to be positioned inside the 2D plane. 3DLandscape is only available with a color acquisition image. 3DLandscape is greyed out in the Menu Select pull-down menu with Black-and-White mode acquisition image.
Render	Changes between the rendered image view and the view of a volume of interest. The volume of interest shows the acquired ultrasound images transformed into an isotropic rectangular coordinate system. The volume of interest can be manipulated as described above.
Reslice	Selections: Cube, Virtual Rescan, and Cubic Plane. Cube. The VOI shows the acquired ultrasound images transformed into an isotropic rectangular coordinate system. This mode allows you to work simultaneously with six cut planes. Virtual Rescan. The marked cut planes under Reslice Cube (red border) is displayed without any perspective distortions, e.g., parallel to the screen. This allows you to move through the volume one slice at a time in any direction. Cubic Plane. Only one cut plane view is shown in a perspective displayed VOI. The cut plane can be moved freely without any limitations.
Overwrite User1/2/3	Overwrite the application presets defined in 3D Acquisition.
Comments	To initiate the comment mode on the 3D image.

Movie 3D

Descriptions and instructions for using Movie 3D follow

Table 6-19 Movie 3D Descriptions and Instructions for Use

3D Parameter	Description
Manual	An animated rotation of the rendered image can be calculated and displayed by this function. Using this function, you first need to define the start and end position of the rotation. To define this, move the VOI to the start position, the press Define Start. Move the VOI to the end position and press Define End.
Movie 360	The calculation and display of a complete rotation around the axis, defined by the Axis button, starts in steps of 15 degrees.
Auto Movie	Initializes the calculation and display of a 3D movie. A rotation of 30 degrees left and right around the actual image position (either the default position after acquisition or the position that was manually defined by manipulating the 3D volume of interest) is shown. For this 60 degree rotation, eleven images in steps of 6 degrees are calculated.
Axis	All rotations (Auto Move and Movie 360) are calculated as rotations around the specified axis (X, Y, or Z).
Movie Speed	You can adjust the speed of any 3D rotation.
Pause	Stops and restarts the rotation. As soon as Pause is pressed, the different rotation steps can be displayed by moving the Trackball.
Colorize/Contrast	Colorizes the 3D render or adds contrast to the 3D rendered image.

Elastography (Option)


Elastography shows the spatial distribution of tissue elasticity properties in a region of interest by estimating the strain before and after tissue distortion caused by external or internal forces. The strain estimation is filtered and scaled to provide a smooth presentation when displayed.

Below is an example of Elastography. The Elastography color map/bar appears, this image appears in dual image mode, and Elastography imaging parameters appear on the right side of the screen as E.

Elastography can be activated via the User defined key on the Control Panel or on the touch panel.

Elasto is only available on Linear and Convex probes.

Figure 6-30 Elastography Example



- 1. Color Map
- 2. Color Bar

Activating Elastography

Table 6-20 Elastography Parameters Description

Parameter	Description
Axial Smoothing	Controls the smoothness of the elastography image in the axial direction. A higher value means a smoother image.
Lateral Smoothing	Controls the smoothness of the elastography image in the lateral direction. A higher value means a smoother image.
Window	Controls the RF data segment size for the motion tracking. A higher value of Window gives a better signal to noise ratio (SNR) at the cost of axial resolution.

Parameter	Description
Map	Controls the elastography maps. Six different maps are available with various contrast and color schemes, including a gray scale map.
Frame Average	Controls the persistence of the elastography images.
Frequency	Controls the transmit frequency.
Soft Compress	Individually controls the image enhancement for the softer than average tissues.
Hard Compress	Individually controls the image enhancement for the harder than average tissues.
Scale	Controls the time interval between consecutive firings. A lower value dictates a higher sensitivity to weak manual motion.
Transparency	High values bring out the tissue behind the elastography data. You adjust via the Color Gain control; this imaging parameter appears as a “T” on the right-hand portion of the display.
Biopsy Kit	Biopsy Kit.
Sample Vol	Controls the transmit pulse length. A higher value means longer transmit pulse which gives better SNR but reduces axial resolution.
Frame Reject	Controls how many frames get rejected due to low quality vertical motion. A higher value means more frames get rejected. A rejected frame has a completely transparent ROI with the B-Mode background showing through.
Noise Reject	Controls how many frames get rejected due to lateral and elevational motion. A higher value means more frames get rejected. A rejected frame has a completely transparent ROI with the B-Mode background showing through.
Line Density	Optimizes B-Mode frame rate or spatial resolution for the best possible image.

How to Use

The Elastography image is achieved by pulsating the probe while you are scanning the anatomy of interest. Here are some criteria to use:

Handheld elasticity imaging can be very dynamic as the size of distortion depends on the movement of the hand-held probe. To maintain stable and consistent displayed strains, pay attention to the Quality graph. Two forms of feedback are provided. In either form, an ideal manual compression is indicated by a high value feedback. In addition, apply the following post-processing controls: Smoothing, Window, Scaling, and Frame Averaging.

Elastography displays harder tissue in blue and softer tissue in red. To enhance Blue area, increase Hard Compress; to enhance Red area, increase Soft Compress. To enhance elastography contrast, reselect the Color Map.

If you need more resolution, increase Smoothing, increase Frequency, reduce Sample Volume, or reduce Window.

If you need a smoother image, increase Window, Axial Smoothing, Lateral Smoothing, and Sample Volume.

If the images seem too flashy, decrease Frame Reject to 1.0 and Noise Reject to have consistent imaging throughout.


Additional Notes for Elastography

NOTE

Elasticity Index is not available in the United States.

Definition. Elasticity index (E) is defined as the value from 0.0 to 6.0. This index shows the color distribution within the measured circle (ROI) relative to the whole ROI box. A higher value means “higher stiffness” and “blue color is dominant”.


Figure 6-31 Elastography Quantification Definitions



Elasto Raw Data has two data sets of index and color for each pixel. Therefore, this index is not affected by the color appearance change (color map and hard/soft compression). Index detects smaller strain difference of hard region by assigning wider dynamic range for a smaller strain.

2D Measurement. Measures the elasticity index of a single ROI and calculates the ratio between two ROIs. You can label two different regions. You can also set the ROI on a reference screen (B mode) with dual screen.

Figure 6-32 2D Measurements



HINTS

- This is a relative quantification tool based on freehand manual palpation technology. It cannot show the stiffness by the kPa (kilopascal).
- There is no compatibility among manufacturers regarding the value. It depends on their strain imaging technology and definition of the value.
- Colors indicate degree of stiffness and do not directly correlate to a specific tissue type. Interpretation of what the tissues are and how to apply these ratios clinically is at the discretion of the user.
- Elastography physics dictates that cystic structures will be displayed with a three-layer pattern. This three-layer pattern will start with blue on the ultrasound system factory default map (which corresponds to hard), then progresses to green and then to red (which corresponds to soft). The posterior displacement of elastography patterns also may cause the B-Mode cyst to consist primarily of blue with the green to red being posterior to the B-Mode cyst. You need to be aware of the three-layer pattern of a cyst in elastography. Utilizing elastography quantification and setting the ROI in the blue portion of the three-layer pattern on the cyst and then setting the ROI in the “normal” tissue may cause you to misinterpret the elastography quantification ratio as the cyst to be hard as compared to the “normal” tissue.

Contrast Imaging (Option)

NOTE

Please be advised that Contrast imaging MAY NOT BE available on your system. Contrast agents for radiology use are not yet available.



CAUTION

Contrast Imaging is not available in USA.



WARNING

Appropriate training

Only physicians or echo technicians who have received appropriate training can use the Contrast applications.



WARNING

Please use appropriate contrast agent.

Misdiagnosis based on image artifacts.

Misdiagnosis in ultrasound contrast images may be caused by several artifacts, most importantly:

Motion artifacts: gives rise to signals independently of contrast presence. This may be caused by patient movement; including respiration, or by probe movement influenced by the operator.

Regional drop outs: caused by unintentional destruction of the contrast agent, too low concentration of contrast agent, poor acoustic penetration due to rib/lung shadows or system failing to detect the contrast agent due to erroneous settings induced by the operator.

Tissue harmonics: gives contrast-like signals independently of the presence of contrast agent.



CAUTION

Cardiac rhythm disturbances during perfusion studies using gas ultrasound contrast agents have been observed in the diagnostic range of Mechanical Index (MI) values and that, for details, to see the specific package insert for the contrast agent being used.



WARNING

Read and follow contrast agent instructions provided by the manufacturer.

NOTE

The ultrasound system is designed for compatibility with commercially available ultrasound contrast agents. Because the availability of these agents is subject to government regulation and approval, product features intended for use with these agents may not be commercially marketed nor made available before the contrast agent is cleared for use. Contrast-related product features are enabled only on systems for delivery to an authorized country or region of use.

By adjusting the acoustic output, you can enhance either contrast harmonics or stimulated acoustic emission (SAE).

Benefits

Injected contrast agents re-emit incident acoustic energy at a harmonic frequency much more efficiently than the surrounding tissue. Blood containing the contrast agent stands out brightly against a dark background of normal tissue.

Clinical Use

Possible clinical uses are to detect and characterize tumors of the liver, kidney, and pancreas and to enhance flow signals in the determination of stenosis or thrombus.

Affect on other controls

The default acoustic output adjusts for contrast imaging and the Power Output key provides more subtle gradations for use while in contrast imaging. When you exit Contrast Imaging, the system returns the acoustic output to its original setting. When you reactivate Contrast Imaging, the system enters the default Contrast Mode.

Most system controls are available (Depth, Zoom, Colorize, etc.). However, some controls are not available (Anatomical M Mode and Rejection).

Controls adjusted while in Contrast Imaging retain these values when you exit Contrast Imaging (except for post-processing controls).

Bioeffect

Activating Contrast Imaging may change the TI and/or MI. Observe the output display for possible effects.

Feature Availability

3D is available; Multi Image and LOGIQView are not available.

Probe

Contrast is available on 4C-RS and C1-5-RS probes in abdomen and pediatric applications, L3-12-RS probe in vascular / Peripheral Vascular and Small Parts applications and IC9-RS probe on gynecology and urology applications.

Mode**Reference Mode**

Reference (Ref) Mode is to image the anatomical reference, not the contrast enhancement.

Contrast Mode

There are several contrast imaging techniques. Note that the appropriate imaging technique may vary by agent and application. In other words, the imaging technique is not dedicated for the agent and vice versa.

Parameters

You can be preset almost controls in the Utility --> Imaging --> CON or Ref Tab.

Single View/Dual View

Description

Single View displays the active Contrast image only.
Dual View displays the active Contrast image and the reference image simultaneously.

Visualization

Description

Define the display technique.

Values

- Contrast. Displays the contrast-enhanced image.
- Tissue. Displays the tissue image.
- Hybrid Contrast. Displays the contrast-enhanced image and the tissue image using Hybrid Map.

Visualization/Map Indication

Visualization/Map Indication

Figure 6-33 Imaging parameter

CON	
Frq	Gen
Gn	45
S/A	0/4
Map	K/0
D	17.0
DR	57
AO%	-
Trig	0-1
Vis	C
Tch	AM

1

2

1. Gray/Colorize Map for current visualization.
2. Visualization: C=Contrast, T=Tissue, HC=Hybrid Contrast

Hybrid Map

Select the hybrid map for the Hybrid Contrast visualization in the Dual/Hybrid Display.

Gray Map

You can select a map for the contrast image and reference tissue image independently.

Adjusting

Press **Gray Map** and select the map in each tab.

Contrast Color

Select the color map for the Contrast mode image.

Adjusting

Press **Colorize** and select the map.

Contrast only

Only transmit and receive the acoustic power signal for the contrast-enhanced image.

Adjusting

Press **Contrast Only**.

Frequency

Select frequency type for each contrast technique.

Adjusting

Frequency can be preset in the *Utility -> Imaging -> CON*.

- Frequency (AM): Res, Gen, Pen

Sonazoid

If you use Sonazoid for Contrast agent, you have to check Sonazoid in ***Utility > Imaging > CON***.

Target MI

Target MI Control provides the automated adjustment of the acoustic output to keep the specified target MI value to reduce the unexpected results on the contrast exam.

Target MI control can be preset in the ***Utility > Imaging > CON and Ref***.

Accumulation

Accumulation enhances the flow in an image.

If Accumulation is turned off, then Frame Averaging is used; if Accumulation value is set, then Accumulation is used.

Availability

Available in Contrast, Color Flow, PDI, and B-Flow.

Benefit

Accumulation detects the maximum signal and holds it for the level specified.

Max (Maximum) Enhance

Sets the acoustic output to its maximum setting (100%)

Values

On/Off. When you deactivate Max Enhance, the acoustic output is returned to its previous setting. Max Enhance is deactivated by the system when you turn it off, change probes, or change the contrast technique.

Benefits

This control provides quick transition to High MI imaging. This allows the user a quick one-button push to destroy the agent. Useful when the user is interested in the bubble wash-in characteristics of the anatomy being scanned.

Trigger

The Contrast Trigger scans images at set intervals, delaying imaging according to the time delay that you specify.

Adjusting

On/Off. Press **Trigger** assigned rotary.

Time Delay

The Contrast Trigger scans images at set intervals, delaying imaging according to the time delay that you specify.

Rotate **Time Delay** assigned rotary.

Adjusting

Set the default seconds of Time Delay in **Utility > Imaging > Con.**

Flash

Description

This feature provides a way to expose the higher acoustic power for a specified time duration by pressing a control once.

Adjusting

Select **Flash** on the Touch Panel.

NOTE

Set the frame numbers to scanned with the higher acoustic power in **Utility > Imaging > Con > Flash Frames**. The frame numbers defined by Flash Frames are applied to both the contrast imaging modes.

- If you select Flash once in the Con menus, the system scans with 100% acoustic output burst pulses for the specified number of frames. The acoustic output then reverts to the original settings.
- When Max Enhance is ON for the contrast imaging modes, the system keeps Max Enhance = ON with no acoustic output change when Flash is selected in the Con menu.

Contrast Clock (Timer)

You can use the Contrast Clock by activating it at the time of injection and deactivating it at the end of the exam.

Two timers can be displayed on the bottom-left corner in image area and info area for several injection.

NOTE

You can also configure the system to perform a countdown for the contrast injection with the **Utility > System > System Imaging > Countdown Time** for Contrast preset.

Adjusting

1. Press **Contrast Clock 1** to start/stop the T1 timer.
2. Press **Contrast Clock 2** to start/stop the T2 timer.

Display

There are two areas on the screen where the Contrast Clock displays: on the image and on the lower, left-hand portion of the display. The timer on the image freezes when you freeze the image (the timer updates when you unfreeze the image). However, the timer located on the lower, left-hand portion of the display continues to display over a freeze, probe change, mode change, multi image, and zoom.

The timer also appears on CINE Loops and archived images.

Benefits

The Contrast Clock measures the time since injection.

You can save the data of contrast clock to an external file by using Export Traces of TIC.

1. Press **Freeze**. Scroll with the Trackball to show Cine tab.
2. Press **TIC Analysis** on the Touch Panel to enter TIC application.
3. Put a ROI on the image.
4. Press Export Traces. Type the file name and store it to the storage device.

Relationship with other controls

- L + R Simultaneous Display (L: Tissue, R: Active Visualization)
 - Accumulation/Cine Capture
 - Applied on the right side image only.
 - You cannot compare On and Off image using L + R.
 - SRI-HD
 - Applied on both side image.
 - Frame Average
 - Applied on both side image.
- TIC
 - Measured on active visualization, except "Hybrid Contrast".
 - "Hybrid Contrast" is disabled and measured on "Contrast" visualization.
- Easy 3D

- Build volume data from active visualization, except “Hybrid Contrast”.
- Hybrid map is disabled and forces “Contrast” visualization.
- Archive
 - The raw data size increases between two to three times the size compared to the previous raw data. It takes a longer period of time to save as Cine.

Time Intensity Curve (TIC) Analysis

The basic TIC process works as follows

1. Scan the patient after injecting the contrast agent.
2. Watch the agent flow through the anatomy of interest.
3. When the desired contrast effect has been visualized, freeze the image and select a range of images for analysis.
4. Position an ROI (region of interest) on one of those images where the contrast effect is visible.
5. The system then calculates the mean pixel intensity within that ROI for all frames in the user designated loop and plots the resulting data as a function of time.

You can also choose to fit this data to one of several mathematical functions. The fundamental idea is that the contrast effect flowing through the organ of interest can be modeled mathematically, and details of the wash in and washout of the agent can be gleaned by analyzing the numerical parameters of the mathematical model.

Activating TIC

1. Scan the patient in Contrast mode or select a desired cine loop from the stored images.

NOTE

Images from the current scan session (already in CINE) or from a saved image loop can be used for TIC analysis.

NOTE

TIC analysis is only available if the user has selected an image loop. If the user has selected a saved still image (just one frame), TIC analysis is not available.

2. Press **Freeze**. Move the Trackball

NOTE

The TIC package is only available when the system is in FREEZE mode.

3. **TIC Analysis** displays on the Cine Touch Panel.

NOTE

Subsequent TIC processing is done on the images stored in CINE memory.

4. Select **TIC Analysis**. The TIC Analysis screen and Touch Panel displays. To toggle the trackball function between QA and Scroll, press **Scan Area** key.

Figure 6-34 TIC Touch Panel




Table 6-21 TIC Touch Panel Description

Parameter	Description
Max Gradient	Between the CINE Start and End frame, displays the time and gradient that becomes the maximum gradient.
Accumulation	Enhances the flow in an image.
Disable Frame	The current frame is excluded from the CINE Loop display.
Enable All Frames	Re-enables disabled frames.
Horizontal Sweep	Allows you to increase/decrease the time interval over which to plot the TIC curve.
Smoothing	Smooths the trace displayed by applying a filter over a defined time window. Both the filter type and time window are user-selectable. The type of filter available depends on the analysis signal displayed.
Delete Sample Area	Removes selected sample area from the CINE Loop window and accompanying trace in the Analysis window. The Trackball marker must be pointed at an anchored sample area.
Export Trace	Saves trace data in ASCII format, readable in spreadsheet programs. If present, trace data for physiological traces are also exported.
Curve Fitting	Toggles between Wash-In, Wash-Out, and Off.
Set Sample Area Shape	Enables resizing of a selected sample area by setting height, width, and tilt angle. The Trackball marker must be pointed at an anchored sample area.

NOTE

The CINE Touch Panel is available to use when in TIC Analysis mode. CINE Touch Panel controls retain their traditional functions. While the CINE Touch Panel is on the top level, the Trackball and other imaging controls perform as they usually do in standard CINE mode.

Exiting TIC Analysis

There are several methods to exit TIC Analysis.

- Select **Exit TIC Analysis** on the TIC Touch Panel.
- Press **Freeze** to unfreeze and resume scanning.

TIC Analysis Screen Description

Figure 6-35 TIC Analysis Screen




Table 6-22 QAnalysis Screen Description

1.	Contrast Cineloop Window Sample Area: Indicates sampling position of the intensity (contrast) trace. The sample area is color-coded: the first sample area is yellow, the second green, etc.
----	---

2.	Displays time-intensity curve. • Y axis: Intensity scale (logarithmic) (db) or linear acoustic units (AU). • X axis: Time(s) or Dt(s), elapsed time from previous frame. • ECG (where available -- not shown), Frame Marker: displays ECG trace (where available), the current frame marker and the start and stop markers for the cineloop. • Time at cursor position and velocity at cursor position. • Intensity (dB or AU) at cursor position. • Intensity (dB or AU) at frame marker position (color coded)
3.	Sample Area Tools. • Pencil Icon: Creates a sample area based on freehand drawing. • Shape Icon: Creates a sample area with a pre-defined circular/ellipse shape.

System Menu

The System menu on the monitor display can be used instead of the Touch Panel.


Analysis Window System Menu

Position the cursor over the analysis window and press the **System Menu** rotary. The Analysis Window system menu displays at the cursor position.

NOTE

The system menu is dependent on mode.


Figure 6-36 Analysis Window System Menu - example



ROI System Menu

Position the cursor on the ROI and press the **System Menu** rotary. The ROI system menu displays at the cursor position.


Figure 6-37 ROI System Menu



Frame Marker System Menu

Position the cursor on the Frame Marker in the Analysis Window and press the **System Menu** rotary. The Frame Marker system menu displays at the cursor position.

Figure 6-38 Frame Marker System Menu



Selecting TIC Analysis Image Range

A range of frames is selected for the TIC analysis in Cine mode (before accessing TIC Analysis). Only the frames in this range are used for the TIC analysis.

If a range is not selected prior to accessing the TIC Analysis, the system uses the default Cine start and end frames as the default TIC start and stop frames.

1. The first frame in the analysis series is selected by adjusting the **Start Frame** control to the desired frame
OR
using the **Trackball** or the **Frame by Frame** control to select the desired first frame and then selecting the **Start Frame** control.
2. The last frame in the analysis series is selected by adjusting the CINE **End Frame** control to the desired frame
OR

using the **Trackball** or the **Frame by Frame** control to select the desired last frame and then selecting the **End Frame** control.

Generating a Trace

Up to eight traces can be generated.

About the sample area

The sample area can be in three different states:

- Free sample area: freely moving sample area (QA cursor) before anchoring.

NOTE

The free sample area disappears when the QA cursor is moved over a static anchored frame.

- Static sample area: the free sample area is anchored by pressing Set.
- Dynamic anchored sample area: the sample area is anchored in two or more frames (see Manual tracking below). In these particular frames, the sample area is displayed with an anchor. The sample area moves smoothly between the anchored positions when playing/scrolling the cineloop.

Trace from a pre-defined sample area

1. Press the **Scan Area** on the control panel to select QA function for trackball.
2. If necessary, select the sample area Ellipse ROI button (shape icon on the monitor display).
3. Move the cursor to one of the Cineloop windows using the **Trackball**.

The cursor is changed to a sample area (white circle). A preview of the trace is displayed in the Analysis window.

4. Press **Set** to anchor the sample area.

In this frame, the sample area is marked with an anchor. If the cineloop has more than one heart cycle, a sample area will also be anchored in the corresponding frame in the next heart cycle.

The trace is updated accordingly in the Analysis window.

Trace from freehand sample area

1. Select the Freehand ROI button (pencil icon on the monitor display).
2. Move the cursor to one of the Cineloop windows using the **Trackball**.
3. Trace the outline of the desired ROI by moving the caliper with the **Trackball**.
4. Press **Set** to anchor the sample area.

The sample area is automatically closed and the trace is updated accordingly in the Analysis window.

Manual tracking of the sample area (dynamic anchored sample area)

The sample area can be moved within the loop to ensure that data in the trace is generated from the same anatomical location during the cyclic motion of the heart.

1. Place a sample area over a region of interest. Note the anatomical location of the sample area.
2. Scroll to a new frame using the **Trackball**.
3. Press the **Scan Area** on the control panel to select QA function for trackball.
4. Move the cursor to the sample area using the **Trackball**.
5. Press **Set**. The sample area is unanchored.
6. Drag the sample area to the corresponding anatomical location in the new frame.

When the sample area is anchored in more than one frame, linear interpolation is performed so that the sample area is smoothly moved between the anchored positions in the selected frames when running the cineloop.

NOTE

In the original frame and this particular frame the sample area is marked with an anchor.

7. Press the **Scan Area** on the control panel to select Scroll function for trackball.
8. Using the **Trackball**, scroll through the cineloop and control that the sample area follows the moving anatomical structure.
9. Add anchored sample areas in several frames to obtain a more accurate displacement of the sample area.

Moving a dynamic anchored sample area

1. Freeze the image.
2. Press the **Scan Area** on the control panel to select Scroll function for trackball.
3. Using the **Trackball**, browse through the cineloop to display one of the frames where the sample area was anchored.

NOTE

In these frames, the sample area is marked with an anchor.

4. Press the **Scan Area** on the control panel to select QA function for trackball.
5. Move the cursor to the sample area using the **Trackball**.
6. Press **Set**. The sample area is unanchored.
7. Drag the sample area to a new location.
8. Press **Set** to anchor the sample area to the new location.

If you want to move the sample area to the same depth, select Move (same depth) from the System Menu.

Zooming in the Analysis window

To zoom:

1. In the Analysis window, press and hold down the **Set** key while dragging the cursor to define the zooming area.
2. Release the **Set** key.

To unzoom:

1. Press the Cursor key in the Analysis window. The system menu displays.
2. Select **Unzoom**.

NOTE

Shown only in zoom mode.

Delete a trace

The user can delete all traces at once or one at a time.

1. If necessary, press the **Scan Area** on the control panel to select QA function for trackball.
2. Move the cursor over one of the sample area. Confirm that cursor is changed to hand icon.
3. Press the **Delete Sample Area** on the touch panel.
4. Select **Current Sample** or **Delete all** as necessary.

Disabling/Enabling the frame

Frame disabling excludes the actual frame from the cineloop display. Frame disabling is available only with contrast data.

Disabling the frame from the frame marker

To disable One Frame :

1. Use the **Trackball** to move the cursor to the frame maker to disable.
2. Press **Set** to disable the frame.
3. The frame marker is changed from green to red to indicate the frame has been disabled.

NOTE

The disabled frame is no longer displayed in the reference window when scrolling through CINE memory.

Disabling multi-frames from the frame marker

1. Use the Trackball to move the cursor to the first frame maker to disable.
2. Press and hold down **Set**
3. Move the cursor with the Trackball to the last frame to be disabled and release Set.

The marker is turn red and the data from that frame is removed from the trace and any subsequent trace processing.

Disabling a frame from the cineloop window

1. Use the Trackball to move the cursor to the cineloop window.
2. Press the **System Menu** rotary. The system menu displays.
3. Select **Disable frame**.

The current frame is disabled and the corresponding frame marker displays red.

Disabling ECG triggered frame (where available)

In a multi-cycle acquisition, the user may deselect all frames in all heart cycles but a selected one. This function can be used for example to select a particular systolic frame for each heart cycle.

1. Scroll through the cineloop to identify the cardiac phase to analyze or identify the cardiac phase on the ECG trace (where available).
2. Position the cursor on the ECG trace (where available) and press **System Menu** rotary. The system menu displays.
3. Select **ECG triggering** (where available).


All frames in all heart cycles are disabled except for the selected and corresponding frames in the other heart cycles.

To enable the frames

To re-enable all deleted frames :

1. Position the cursor on the Frame Marker line and press the **System Menu** rotary. The system menu is displayed at the cursor position.
2. Select **Enable all frames**.
3. All disabled frames are re-enabled.

Figure 6-39 Frame markers



4. Analysis Window
5. Frame markers axis
6. Enabled frame (Green)
7. Disabled frame (Red)
8. ECG (where available)
9. Current frame

Manipulating the Sample Area

Up to eight ROIs can be saved on the reference image, with the corresponding eight traces plotted simultaneously on the graph. Each ROI display has a different color, and its corresponding trace data is plotted using that same color.

Once eight ROIs have been saved, the system does not automatically generate an active ROI when the cursor is positioned over the displayed reference image.

The saved ROIs can be a mixture of elliptical and freehand ROIs.

When the user repositions an ROI, the old trace data is erased from the plot and the trace data for the new position replotted.


If the ROI position on the last frame of the selected image range is moved, the corresponding ROIs on all frames are repositioned to match the last frame.

The user shall also have the capability of setting separate ROI positions on different frames of the contrast images, and the system shall linearly interpolate the ROI positions for the frames in between the selected frames.

Setting the default sample area shape

1. Select **Set sample area shape**. The Information Box displays.

Figure 6-40 Sample Area Information Box



2. Select Height, Width and Tilt angle.
3. Select **Set as default**. The current ROI size is set as the default for subsequent Ellipse ROIs.

Reshaping a Sample Area

To reshape the sample area:

1. Position the cursor on the ROI to reshape and press the **System Menu** rotary.
2. The ROI system menu displays. Select **Set sample area shape**.
3. Adjust Height, Width and Tilt angle.
4. Press **OK**. The selected ROI size changes.

Labeling a Sample Area

The sample area label is used to identify data associated with the sample area when exporting.

1. Position the cursor on the ROI to label and press the **System Menu** rotary.
2. The ROI system menu displays. Select **Label sample area**. The Label Dialog box displays.
3. Enter a name for the sample area.
4. Select **OK**.

Sample Area Shapes

There are two different methods for determining the shapes of the sample area.

Ellipse ROI

1. Select the ellipse icon (shape icon on the monitor display).
2. When the trackball positions the image display cursor over the reference image(s), an elliptical ROI is automatically generated and displays on the reference image(s).
3. The average intensity value inside the ellipse is calculated for every image in the image analysis range and plotted in the image display area.
4. The last generated or selected ellipse is considered the active ROI, and its trace plot automatically updates as the user repositions it on the reference image. Old traces are erased.
5. When scanning with an elliptical ROI, press **Set** to fix the ROI position and freeze its corresponding trace on the plot. A new active ROI is generated whose position is manipulated by the trackball and whose time-intensity curve traces will be plotted as before, while the previous ROI and trace remain fixed at the points they were saved at.

NOTE

Elliptical ROIs can be positioned in any manner that keeps their center within the image boundaries. In the case that part of the ROI is outside the image boundary, only data from within the image boundary is used for calculating the mean intensity value.

NOTE

You can change the size of the Ellipse ROI by adjusting the Ellipse control.

Freehand ROI

1. Select Freehand icon (pencil icon on the monitor display).
Use the **Trackball** to position the caliper on the reference image at the start point. Press **Set** to fix the start point.
2. Trace the outline of the desired ROI by moving the caliper with the **Trackball**.
3. When a suitable ROI has been drawn, press **Set** a second time.

The system automatically links the start point to the end point by drawing a straight line between them. The caliper is then free for repositioning for another freehand ROI.

NOTE

You cannot go outside the image boundary when drawing a freehand ROI.

Copy, move and paste a Sample Area

To copy and paste the ROI,

1. Move the cursor over the ROI and press the **System Menu** rotary. The system menu displays.
2. Select **Copy sample area**.
3. Move the cursor to the desired location for the copied ROI and press the **System Menu** Key. The system menu displays.

4. Select Paste sample area.

To copy and move the ROI,

1. Move the cursor over the ROI and press the **System Menu** rotary. The system menu displays.
2. Select **Copy & move**. Or if you want to move to the same depth as the original ROI, select **Copy & move (same depth)**.
3. Move the copied ROI using the **Trackball**. Press **Set** to fix the position.

Deleting a Sample Area

Sample ROIs and their corresponding traces can be deleted using **Delete Sample Area**.

1. Select **Delete Sample Area**; a pull-down menu displays.
 2. Select **Current sample** to delete the currently active ROI.
- Select **Delete all** to delete all currently set ROIs and all of their traces.

NOTE

The corresponding traces for the deleted ROIs are erased from the plot.

NOTE

Deleting an ROI causes the ROIs to be deleted from all frames in the analysis loop.

TIC Plot Control

The following controls are user configurable presets which are configurable via the Utility Menu or through the pull-down menu in TIC Analysis mode. When using the pull-down menu:

1. Position the cursor over the analysis window and press the **System Menu** rotary. The system menu displays at the cursor position.
2. Select the appropriate parameter.

Vertical Unit


When analyzing the contrast data, the Y-axis can be set to display either logarithmic scale (dB) or linear, acoustic units (AU) for both tissue intensity (2D) or Angio intensity data.

To toggle between dB and acoustical display units for the Y-axis.

- **dB**—The traditional log compressed B-mode data is used to calculate the time-intensity curve values.

- **Acoustic**—The system reverse the log compression function to provide un-log compressed data for the TIC analysis.

Figure 6-41 Vertical Unit Pop-up menu




Vertical auto-scaling

The system can be configured to display the full unit range or a range according to the maximum and minimum values of the displayed trace(s) (auto-scaling function). In addition, the auto-scaling function can be set to be live update (updates while the sample area is moved) or delayed (updated when the sample area is anchored).

- **Delayed**—The system automatically rescale the vertical axis of the trace graph only when a new ROI is saved, to account for changing input dynamic range.
- **On**—The system automatically rescale the vertical axis of the trace graph every time the currently selected (active) ROI is moved.
- **Off**—Disable any automatic scaling of the vertical axis. There is user-defined system defaults on the system preset page for the fixed vertical scale to be used for the plot.


Figure 6-42 Vertical Autoscale Pop-up menu



Line Style

- **Solid**—Setting the results in a plotted trace that does not display small boxes at the data points
- **Squares**—Setting the results in a plot where small squares are displayed at each data point, and the squares are linked together by lines.

Figure 6-43 Line style Pop-up menu



Horizontal Sweep

Horizontal Sweep allows you to increase or decrease the time interval over which to plot the TIC curve.

The default is the user selected image range. If the user has not yet selected a first and last frame, the first and last default frames from the displayed CINE loop are used.

Smoothing

The system can smooth the traces displayed by applying a filter over a defined time window. The type of filter available is depending on the analysis signal displayed.

1. Select **Smoothing**.

OR

Position the cursor over the analysis window and press the **System Menu** rotary. The System menu is displayed at the cursor position. Select **Smoothing**.

NOTE

When smoothing is turned on, it applies to all traces in the plot window.

2. The smoothing filter list displays. Select the appropriate parameter.


NOTE

If 'dt' is selected for Horizontal Scale, you cannot use Smoothing.

Trace Measurement Gradient

Gradient is displayed on the screen instead of Intensity (dB or AU). The gradient calculates from 7 points (includes previous and next frames).


Figure 6-44 Gradient



Show graph (Gradient plot)

1. Gradient Plot displays on the TIC system menu.


Figure 6-45 System Menu



2. Select the parameter.
 - Off: A graph plots TIC.
 - Gradient: Two graphs plot TIC and TIC gradient.
 - Unit of Y-axis is dB or AU in case of intensity.
 - The unit is $d(\text{dB})/dt$ or $d(\text{AU})/dt$ in case of the intensity gradient.


- Gradient values for the current frame are displayed in the upper right corner of the graph.

Figure 6-46 Gradient



- Gradient Derivative: Two graphs plot TIC and TIC gradient derivative.
 - The Y-axis units is $d^2(\text{dB})/dt^2$ or $d^2(\text{AU})/dts$ in case of the intensity gradient derivative.
 - Gradient derivative values for the current frame are displayed in the upper right corner of the graph.


Figure 6-47 Gradient derivative



Advanced Features

- All: Three graphs plot TIC, TIC gradient and TIC gradient derivative.


Figure 6-48 All Gradient



Max Gradient

Displays the time and gradient that becomes the maximum gradient between the CINE start and end frame.

Figure 6-49 Max Gradient



Curve Fit


Curve fitting analysis: for research studies of perfusion rates using contrast agents.

1. Select **Curve Fit**.
OR

Position the cursor over the analysis window and press the **System Menu** key. The system menu displays at the cursor position. Select **Curve Fit**.

2. The Curve Fit selection list displays.


Figure 6-50 Curve Fit Selection List



- **Off**—Remove the fitted curves from the plot and the fit parameters from the display.
- **Wash-in**—Used to find and estimate the local perfusion rate using the contrast agent. Exponential wash-in is described by the function:
 $Y(t) = A(1 - \exp(-kt)) + B$, where:
 - A (dB or AU) is the intensity from the contrast agent.
 - B (dB or AU) is the intensity at time $t=0$ (defined as the time of the left marker). This corresponds to the tissue (baseline) signal if no contrast is present at the selected starting point.
 - $A + B = \text{contrast} + \text{tissue} = \text{plateau level}$.
 - k (1/s) is a time constant.
- **Wash-out**—Used to find and estimate a local wash-out rate. Exponential wash-out is described by the function:
 $Y(t) = A\exp(-kt) + B$, where :
 - A (dB or AU) is the intensity from the contrast agent.
 - B (dB or AU) is the intensity from the tissue = baseline signal.
 - $A + B$ is the initial intensity level.
 - k (1/s) is a time constant.
- **Gamma variate**
 $Y(t) = At^C\exp(-kt) + B$

Parameters of Gamma curve fitting

Figure 6-51 Gamma Curve




- tc : Increasing function ($C > 0$) for “Wash-in”.
For larger C , the intensity increases quickly before the peak.
- $\exp(-kt)$: Decreasing function ($k > 0$) for “Wash-out”.
For larger k , the intensity decreases quickly after the peak.
- B : Intercept intensity at $t=0$.
- The peak intensity of the curve is affected by all parameters.
Larger A , larger B , larger C , and smaller k make larger peak. The peak time is calculated by C/k .
- MSE: Mean Square Error
If the MSE is small, the difference of actual data and the fitted curve is small.

Set Start/End Frame for Curve Fit per ROI position

1. Generate TIC and perform a Curve Fit.

In this state, the Curve Fit graph is drawn from Cine Start Frame to Cine End Frame for all the ROIs.

Figure 6-52 Curve Fit screen




a Cine Start Frame

b Cine End Frame

2. Select the Start Frame as with Cine or use the trackball to select Start Frame and press **Set** with setting Scroll function to trackball by toggle **Scan Area**.
3. Set QA function to trackball by toggle **Scan Area**, move the arrow pointer on the ROI (hand cursor appears) and press **System Menu**. The system menu appears.
4. Select **Set curve fit start frame** from the menu.
5. Select the End Frame as with Cine or use the trackball to select End Frame and press **Set** with setting Scroll function to trackball by toggle **Scan Area**.
6. Set QA function to trackball by toggle **Scan Area**, move the arrow pointer on the ROI (hand cursor appears) and press **System Menu**. The system menu appears.

7. Select **Set curve fit end frame** from the menu. The ROI colored line displays.

Figure 6-53 New Start Frame and End Frame (Example)



a Start Frame

b End Frame

8. Repeat the above procedures as necessary.

The system retains the start/end frame per ROI while TIC is active. Once the TIC menu is closed, the settings are lost.

Display/Hide Calculation Values

You can specify Time to Peak, Area under the curve, Curve gradient, and Arrival time values for Wash-in, Wash-out, and Gamma Variate curve fitting values.

1. Select **Curve Fitting Parameters** on the Touch Panel. The Curve Fitting Parameters dialog appears.

Figure 6-54 Curve Fitting Parameters Dialog

You can select up to 5 parameters in each group.

Wash In	Wash Out	Gamma Variate
☒ A	☒ A	☒ A
☒ B	☒ B	☒ B
☒ K	☒ K	☒ K
☒ MSE	☒ MSE	☒ C
☐ Time To Peak	☐ Time To Peak	☒ MSE
☐ Area Under The Curve	☐ Area Under The Curve	☐ Time To Peak
☐ Curve Gradient	☐ Curve Gradient	☐ Area Under The Curve
☐ Arrival Time	☐ Arrival Time	☐ Curve Gradient
		☐ Arrival Time

Set As Default Save Cancel

2. Select a maximum of 5 parameters to display for each curve fitting.


- Save as default: saves as a system preset.
- Save: saves as temporary.
- Cancel

NOTE

If you select more than 5 parameters and select “Save” or “Save as default,” you will be prompted the message that “More than 5 parameters selected”.

- The selected parameter displays below the image with the curve fitting active.

Figure 6-55 TIC image with all selected parameters



Exporting Traces (Saving the Trace Data)

You can save the trace data to an external file.

- Select **Export Traces** on the Touch Panel to save the trace data.
OR

Position the cursor to the Analysis window and press the **System Menu** rotary. The system menu displays. Select **Export traces**.

- The following window displays.

Figure 6-56 Export Trace Window

- Location: Select Location which to save.
- Filename: Enter the filename. (Only Text)

3. Select **OK** to save the data and return to the TIC Analysis screen.
 - All displayed ROI traces are saved in the exported file.
 - The fit parameters are included in the trace file if the user has done a curve fit.

All plot data (intensity, gradient and gradient derivative) are exported to a text file by “Export Trace”.

Table 6-23 Example of exported file

Time(s):	Trace 1:	Trace 1 dGrad.:	Trace 1 dGrad.
0.00000	-3.97995e+000	-2.15924e+001	8.05159e+001
0.03121	-5.14631e+000	-1.64719e+001	1.74256e+001
0.06242	-5.75798e+000	-1.27675e+001	-7.78004e+001
0.09362	-6.02222e+000	-1.27675e+001	-1.93426e+002
0.12483	-6.11224e+000	-1.44515e+001	-4.17252e+002

NOTE

The Smoothed trace is the one saved if the user has applied a smoothing filter.

NOTE

Only data from the user selected image range is included in the exported trace file.

NOTE

Data for disabled frames are not be included in the exported trace file.

NOTE

No trace results are saved in the standard image database.

NOTE

Trace results are not shown on the Worksheet.

Annotating the TIC Data

The user can annotate both the reference image and the trace plot displays. Use **Comment** key to type the annotation. See Chapter 6 for reference.

Printing TIC Data

Press the appropriate print key in the TIC mode.


The system capture a single still frame which consists of the plot, the reference image and user annotation.

Needle Recognition (Option)

Needle recognition functionality is available with linear probes (12L-RS, L6-12-RS, 9L-RS, L3-12-RS, L4-20t-RS), convex probes (4C-RS, C1-5-RS) and bi-plane probe (E7C8L-RS).

Below illustration **b** illustrates the position between the needle and probe. Needle recognition can only enhance in-plane needles. The needle angle is defined as the angle between the needle and probe surface.

Figure 6-57 Enhanced area and Needle Angle



NOTE
Needle Recognition is only available in B/CF/PDI.

NOTE
The parameters for adjusting can be set in **Utility > System > B**. To set the desired parameters in the Touch Panel for each button.

Adjusting

- Activating Needle Recognition:
Press user defined **Needle** key on the control panel to activate Needle Recognition.
- Adjusting steer direction:

Needle key cycles through Off/Left or Off/Right depending on the setting in **Utility > Imaging > B** tab. If the system is set for Off/Left/Right, press **Needle** key again after activating the Needle Recognition in order to change the steering angle from left to right.

- Adjusting the needle recognition beam angle:

Adjust **Beam Angle** by rotating the Touch Panel control. Max available angles are up to 50 degrees on 12L-RS, 9L-RS, L3-12-RS and L4-20t-RS, 45 degrees on 4C-RS, C1-5-RS and 30 degrees on L6-12-RS, 20 degrees on E7C8L-RS.

Adjust the needle recognition steering angle to form the needle and the beam angle as perpendicular as possible to get the best needle enhancement.

- Adjusting Needle Gain:

Adjust **Needle Gain** by rotating the Touch Panel control. Available Needle Gain ranges from 0 to 100 in the increment of 5.

NOTE

Increasing the Needle Gain enhances the needle visibility while decreasing the Needle Gain decreases artifacts and noise. Adjust the Needle Gain to balance the needle enhancement versus level of artifacts/noise.

- Adjusting Needle Thickness:

Adjust **Needle Thickness** by pressing the Touch Panel control to change the appearance of needle thickness between thick and thin.


NOTE

"Thick" option displays thicker and more prominent appearance of needle while "Thin" option displays thinner but more realistic needle appearance.

NOTE

The beam divergence and effects of beam steering with a transducer may prevent a segment of the needle shaft from showing in the image if the needle angle is too steep.

Figure 6-58 Example of needle that is too steep



- 1 Needle shaft
- 2 Un-enhanced needle shaft
- 3 Needle tip

NOTE

For best results please insert the needle at perpendicular to the dotted guide line (Figure 13-83), or at a needle angle that is slightly less than the Beam Angle.

Figure 6-59 Example of needle that is perpendicular to the guide line





Figure 6-60 Example of needle angle that is slightly less than the beam angle





HINTS

To help verify the location and trajectory of the needle tip, please use needle movement and / or fluid injection.

Make sure the needle is always in the ultrasound plane by slightly moving or tilting the probe to get the best needle enhancement during the needle procedure.

Switch needle On/Off can help to identify artifacts and other not interesting structures.

Needle Recognition values (including enable/disable, Needle Direction, Beam Angle and Needle Gain, etc.) are returned to factory or user preset value when you change: Exam Category, Exam Calcs, or New Patient.

Preset

Needle Directions - preset via **Utility > Imaging > B** Tab.

Beam Angle - preset via **Utility > Imaging > B** Tab.

Needle Gain - preset via **Utility > Imaging > B** Tab.

Needle Thickness - preset via **Utility > Imaging > B** Tab.

Interactive Needle line - preset via **Utility > Imaging > B** Tab.

Benefits

Provides better biopsy needle visualization than normal B-Mode with steered beam angle and post processing.

Bioeffects

Activating angle Needle Recognition may change the TI and/or MI. Observe the output display for possible effect.

Breast Care

Breast Care provides a protocol of breast scanning which helps user to scan breast by 4 Quadrants or 12 O'clock and axillary fossa consistently. User can use 'P1' key to follow the scanning steps that allows user to focus on performing the exam rather than on controlling the system and can help user to increase consistency while reducing keystrokes. The system automatically invokes the correct modes, advances to the next step in the exam. After finishing whole scanning steps, report will be generated.

Breast Care does not perform automatic breast lesion detection or classification. The user is required to detect and diagnose any breast tissue abnormality, as well as verify any outputs of Breast Scanning, Breast Lesion M&A, BI-RADS® classification.

Breast care is divided into 3 steps: prepare the patient information, scanning and report.


Prepare the patient information and scanning

Before Perform the scan, user can set up modes, PW, Elasto which can be included in scanning steps and 4 Quadrants or 12 O'clock, Review Previous Report options in **Utility > System > General > Breast Care**.

Example of 4 quadrants scanning steps:

1. User selects the appropriate transducer for breast imaging and select breast preset.
2. Configure **Breast Care** as a User Configurable key. Press the configured key to start.

Figure 6-61 Breast Care User Configurable key



3. The Breast Care protocol appears on the left side of the display.
Breast scanning is divided into 4 areas: Right Breast, Right Axillary Fossa, Left Breast, Left Axillary Fossa.
4. Users have the choice to follow the sequential evaluation of the respective anatomy, and instructions for which trackball controls are displayed if the user selects to manually categorize a segment as negative or positive for findings.

Advanced Features

The system will visually display purple for a segment if negative selected by user. Positive segment will display yellow if selected by user.

User can manually mark the lesion in the positive area with bright circles **b** via the button **a** refer to step 5.

If a segment is marked positive by the user, the system will return to that segment for further measurement and color Doppler, refer to steps 6-12.

Figure 6-62 Breast scan screen1





Figure 6-63 Lesion marker



5. Press **P1** key, the screen is split in to dual screen, the frozen longitudinal plane scan image is placed on the left side.
6. The live image is on the right side and the user performs transverse plane scan.


Press **P1** key, perform measurements.

Figure 6-64 Breast scan screen 2



The pie chart is used to display the breast scan area and result:

Figure 6-65 Pie Chart in Breast Care



- a. The current scanned area: the outline is highlighted.
- b. Dark blue area: Negative area, there is no finding/lesion with breast anatomy.
- c. Yellow area: Positive area, there is a finding/lesion with breast anatomy, probable pathology.
- d. Nipple

Breast Productivity Package

There are two features related to Breast:

- Breast Measure Assistant - this contains the auto contour feature and auto Height and Length. It also has measurements related to Breast (distance to nipple, ratio)
- Breast Productivity - this includes lesion measurement folders, summary, etc.

This section covers both the Breast M&A Package and Breast Measure Assistant (Auto Contour) features.

Breast Lesion M&A

Breast Lesion M&A allows you to document up to 30 breast lesions for each breast. Lesion Height/Width/Length, Distance to Nipple, and A/B Ratio are available. Distance to Nipple allows you to enter the value. This is not a calculated measurement.

The Measure Assistant Breast (Auto Contour) feature can also be used to automatically detect and outline the breast lesion.

Worksheet and Summary Worksheets show all the documented right/left breast lesions.

From the Small Parts Model, select the Breast Application. Next, select the Right/Left Lesion (Select Rt Side/Lt Side below the Touch Panel).

Table 6-24 Breast Lesion M&A Touch Panel

Preset Parameter	Description
Position	Specify the position of the lesion: Clock position 1-12 O'clock, Areolar, SubAreolar, Axillary, or "-" (default).
Segment	Specify A, B, C, None, or "-" (default).
Return	Press to return to the previous Touch Panel.
Lesion #	Indicates which lesion you are viewing (Lesion # of Total Number of Lesions). Press the left/right arrow to move from lesion to lesion.
L	Lesion Length
H	Lesion Height
W	Lesion Width
Distance to Nipple	Used to manually enter the distance the lesion is from the nipple.
Auto Contour (HxL)	Press to activate the Auto Contour feature, using the height and length.
Auto Contour (HxW)	Press to activate the Auto Contour feature, using the height and width.
Rt or Lt A/B Ratio	Right or Left Lesion A/B Ratio, measured by Area or Diameter.
Composition	Specify the composition of the lesion: None (-), Solid, Cystic, or Complex.
Delete Lesion	Press to delete this lesion.

BI-RADS (ACR)

BI-RADS® (Breast Imaging-Reporting and Data System) is a risk assessment and quality assurance tool developed by American College of Radiology that provides a widely accepted lexicon and reporting schema for imaging of the breast.

Breast imaging studies are assigned one of seven assessment categories from ACR BI-RADS® Atlas 5th edition quick reference ultrasound section (<https://www.acr.org/Clinical-Resources/Reporting-and-Data-Systems/Bi-Rads#Ultrasound>) by Mendelson EB, Bohm-Velez M, Berg WA, et al. ACR BI-RADS® Ultrasound. In: ACR BI-RADS® Atlas, Breast Imaging Reporting and Data System. Reston, VA, American College of Radiology; 2013.

Copyright Notice/Permission

Copyright © 1992, 1993, 1995, 1998, 2003, 2013 American College of Radiology. All rights reserved.

No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopy, recording, or any information storage and retrieval system, without the prior express written consent of the American College of Radiology.

The ACR logo is a registered trademark and service mark of the American College of Radiology. BI-RADS® is a registered trademark of the American College of Radiology.

Disclaimer

The American College of Radiology (ACR) Breast Imaging Reporting and Data System Atlas (BI-RADS® Atlas) is the product of a collaborative effort among members of various committees of the American College of Radiology. BI-RADS® suggests a standardized method for breast imaging reporting. This system is not meant to dictate individual case management decisions.

All referring physicians and radiologists should be aware of the benefits and limitations of the application of imaging techniques. Imaging techniques classified as investigational by the FDA have not been considered in developing the BI-RADS® Atlas; however, study of new equipment and applications should be encouraged. The ultimate decision regarding breast imaging procedures and treatment must be made by the referring physician and radiologist in light of all the circumstances presented in an individual examination.

Warranty

Warranty limitations:

1. ACR BI-RADS® is provided on an “as is” basis, without warranties of any kind, either express or implied, including, but not limited to, warranties of title, or implied warranties of merchantability or fitness for a particular purpose;
2. The use of BI-RADS® is at your own risk;
3. Access to ACR BI-RADS® may be interrupted and may not be error free;
4. Neither the American College of Radiology nor anyone else involved in creating, producing or delivering ACR BI-RADS® contained therein, shall be liable for any direct or indirect, incidental, special, consequential or punitive damages arising out of your use or inability to use ACR BI-RADS®.

Advanced Features

To activate BI-RADS® on the system:

1. When the operator would like to assess the lesion in breast, frozen image or recall a still image with the lesion.
2. Press Measure and selects **Rt Lesion** or **Lt Lesion** in **Measure** Menu.
3. The operator measures L*H*W* manually or select Auto Contour(HxL)/(HxW). (L: length; H: Height; W: width)

NOTE

BI-RADS assessment can only be performed after completion of lesion measurement.

4. Selects BI-RADS® (ACR) on Touch Panel. Then BI-RADS® (ACR) assessment form appears. The operator selects the assessment descriptor in the assessment form.
5. Enter the characterization for the assessment description.

Figure 6-66 BI-RADS® (ACR) Menu

The screenshot shows a touch panel interface for BI-RADS (ACR) assessment. The title bar at the top reads 'Rt Lesion -- BI-RADS(ACR)' with an information icon and a close button. The interface is divided into several sections:

- Masses**: A section with a collapse icon, containing five rows of dropdown menus for 'Shape', 'Orientation', 'Margin', 'Echo pattern', and 'Posterior features'.
- Calcifications**: A section with a dropdown menu for 'Calcifications'.
- Associated Features**: A section with a collapse icon, containing four rows of dropdown menus for 'Associated features*', 'Skin changes', 'Vascularity', and 'Elasticity assessment'.
- Assessment**: A section with a row of radio buttons for assessment levels: '-', '0', '1', '2', '3', '4A', '4B', '4C', '5', and '6'. The '0' button is currently selected.
- Comments**: A text input field at the bottom for entering additional information.

NOTE

At any time, the operator can click the icon at the right corner of BI-RADS (ACR) form, the information window from the BI-RADS (ACR) appears. To close the BI-RADS (ACR) information window, the operator clicks the icon again.

Figure 6-67 Summarized User Prompt Window

The screenshot shows the BI-RADS (ACR) form interface. On the left, a sidebar titled 'Rt Lesion -- BI-RADS(ACR)' contains expandable sections: Masses, Calcifications, Associated Features, and Assessment. A blue box highlights an information icon (i) in the top right corner of this sidebar. The main area displays the 'ACR BI-RADS® Atlas Fifth Edition-QUICK REFERENCE-Ultrasound' form. This form includes sections for Masses, Calcifications, Associated Features, and Assessment. The Assessment section shows a row of buttons for categories 0 through 6, with category 0 currently selected. The form also includes a 'Comments' field at the bottom.

6. The operator selects BI-RADS® assessment score.

NOTE

The determine definitive diagnosis depends on FNA (fine needle aspiration).

7. Repeat steps 1- 6 to complete multiple descriptor assessments.
8. After the completion of the assessment, operator can press the Report key to review all the measurements and BI-RADS® assessment in report and worksheet.

Table 6-25 BI-RADS® Assessment Categories

Category	Assessment
Category 0	Incomplete - Need additional imaging evaluation
Category 1	Negative
Category 2	Benign
Category 3	Probably benign
Category 4	Suspicious: - Category 4A: low suspicion for malignancy - Category 4B: moderate suspicion for malignancy - Category 4C: high suspicion for malignancy
Category 5	Highly suggestive of malignancy
Category 6	Known biopsy-proven malignancy

9. The operator can select Exit to exit to worksheet and back to scanning page to do scanning or more lesions' BI-RADS® assessment.

Worksheet and Summary Worksheets

Worksheets and Summary Worksheets are provided for all documented Breast Lesions.

To move to the next page, rotate the **Page Change** rotary beneath the Touch Panel.

NOTE

Only defined features are displayed on the Summary Report. To display the undefined features, select **Show Undefined Features** at the bottom of the Summary Worksheet.

Breast Measure Assistant (Auto Contour)

You can request that the system trace/outline the border of a breast lesion using Measure Assistant Breast (Auto Contour). You do this by setting the Region of Interest (ROI) around the lesion; the system can then measure the lesion by drawing the contour around it.

The system can store up to 30 breast lesions. The system tracks these by numbering the lesions consecutively.

To automatically detect the breast lesion on the display,

1. Press Measure.
2. Press Auto Contour (HxW) on the Touch Panel.
3. Place the Cursor in the center of the lesion and press Set. Size the ROI around the lesion. Use the Trackball to resize the ROI.
 - To increase the size of the circle, move the Trackball down and to the right.
 - To decrease the size of the circle, move the Trackball up and to the left.

NOTE

Include the entire lesion, even if additional surrounding tissue is included.

4. Press Set on the Trackball. A trace appears around the lesion.
5. Size the trace via the Trackball.
6. Press Set. The contour around the breast lesion is generated.

NOTE

Multiple breast lesion traces may be generated by the system.

7. Inspect the generated contour for accuracy. If edits are necessary, execute steps 8-10 to edit the contour prior to accepting the measurement. Otherwise, skip to step 11.
8. To cycle through the generated contours, use the Select Contours rotary on the Touch Panel.
9. To edit the selected contour, move the Trackball to appropriately size the edit region and then press **Set** on the Trackball.
10. The blue portion of the contour can be edited by moving the Trackball to the portion of the contour you want to edit.

NOTE

The Caliper closest to the cursor enables editing.

NOTE

To limit the horizontal/vertical editing capabilities, you can set a preset via **Utility > Measure > Advanced > Small Parts > Restrict Breast Contour Caliper Edit**.

11. After you have completed your edits, press **Done** on the bottom Set Key or press **Print** to accept the measurement.

Thyroid Productivity Package

A Thyroid Productivity Package is available.

The Thyroid Productivity Package is only available when the Thyroid Productivity Option is activated.

Table 6-26 Thyroid/Parathyroid/Lymph Node Touch Panel

Preset Parameter	Description
Side	Specify the side: Right, Left, Isthmus.
Worksheet/Summary	Select to view the Worksheet/Summary Worksheet.
Add#1, Add#2, etc.	Cycles through the available lesions, or adds a new lesion/node/nodule, etc.
Rt/Lt Thyroid Rt/Lt Parathyroid Rt/Lt Lymph Node Rt/Lt Nodule	When measuring the Left/Right Thyroid/Parathyroid/Lymph Node/Nodule, these folders highlight on the Touch Panel. Length, Height, and Width are available for all thyroid measurements. The Cortical Thickness measurement is available for the Lymph Node. Show Features is available for all thyroid measurements.
Location	Parathyroid: Specify Upper Gland or Lower Gland Lymph Node: Supraclavicular fossa, Lower cervical, Middle cervical, Upper cervical, Parotid, Submandibular, Submental, Posterior triangle Nodule -- Use Up/Down Toggle to adjust. - Location A: Upper, Lower, Mid, None - Location B: Lateral, Medial, Mid, None
Isthmus Lymph Node Isthmus Nodule	When measuring the Isthmus Lymph Node/Nodule, these folders highlight on the Touch Panel. Length, Height, and Width are available for all isthmus measurements. The Cortical Thickness measurement is available for the Lymph Node. Show Features is available all isthmus measurements.
Show Features - Overall Thyroid	Press to activate the Show Features notations. To add notations for each feature, position the trackball to the right of each feature and press Set. This brings up the available notations. Move the Trackball to highlight a notation and press Set to select a notation. The notation will then appear next to the feature and on the Summary Worksheet. Below is a list of each Feature with its possible notations by measurement type: - Overall Thyroid (Top Level Touch Panel) - Resected: Totally, Partially, None (-) - Appearance: Within normal limits, Abnormal, Symmetric, Asymmetric R>L, Asymmetric L>R, None (-) - Comment

Preset Parameter	Description
Show Features - Lt/Rt Thyroid	Press to activate the Show Features notations. To add notations for each feature, position the trackball to the right of each feature and press Set. This brings up the available notations. Move the Trackball to highlight a notation and press Set to select a notation. The notation will then appear next to the feature and on the Summary Worksheet. Below is a list of each Feature with its possible notations by measurement type: • Lt/Rt Thyroid Resected: Totally, Partially, None (-) Echogenicity: Homogeneous; Coarse; Heterogeneous; Hashimoto, Classic; Hashimoto, Probable; None (-) Vascularity: Normal, Increased, Decreased, None (-) Size: Normal, Enlarged, Small, None (-) Comment Isthmus Comment

Preset Parameter	Description
Show Features - Lt/Rt Parathyroid / Lymph Node / Nodule	Press to activate the Show Features notations. To add notations for each feature, position the trackball to the right of each feature and press Set. This brings up the available notations. Move the Trackball to highlight a notation and press Set to select a notation. The notation will then appear next to the feature and on the Summary Worksheet. Below is a list of each Feature with its possible notations by measurement type: • Lt/Rt Parathyroid -- Upper/Lower Gland Visibility: Visualized, Not Visualized, None (-) Comment • Lt/Rt Lymph Node Appearance: Within normal limits, Suspicious, Pathologic, None (-) Composition: Cystic, Complex, Solid, None (-) Vascularity: Normal, Increased hilar, Increased non-hilar, None (-) Comment • Lt/Rt Nodule Shape: Round, Oval, Irregular, Lobulated, None (-) Margin: Well Defined, Well Defined with halo, Well Defined with partial halo, Well Defined with complete halo, Poorly defined, Irregular, None (-) Composition: Solid, Cystic, Mixed, Complex, Heterogeneous, None (-) Vascularity: Normal, Increased, Decreased, Central Vasc (Avascular, Hypovascular, Isovascular, Hypervascular, Severely Hypervascular), None (-) Calcification: No Calcification, Coarse central, Coarse rim, Punctate scattered, Punctate clumped, Colloid, Mixed, None(-) Comment

Preset Parameter	Description
Show Features - Isthmus	Press to activate the Show Features notations. To add notations for each feature, position the trackball to the right of each feature and press Set. This brings up the available notations. Move the Trackball to highlight a notation and press Set to select a notation. The notation will then appear next to the feature. Below is a list of each Feature with its possible notations by measurement type: - Isthmus Lymph Node Appearance: Within normal limits, Pathologic, None (-) Composition: Cystic, Complex, Solid, None (-) Vascularity: Normal, Increased, None (-) Comment - Isthmus Nodule Shape: Round, Oval, Irregular, Lobulated, None (-) Margin: Well Defined, Well Defined with halo, Well Defined with partial halo, Poorly defined, Irregular, None (-) Composition: Solid, Cystic, Mixed, Complex, Heterogeneous, None (-) Vascularity: Normal, Increased, Central Vasc (Avascular, Hypovascular, Isovascular, Hypervascular, Severely Hypervascular), None (-) Calcification: Coarse central, Coarse rim, Punctuate scattered, Punctuate clumped, Colloid, Mixed, None (-) Comment
Return	Press to return to the previous Touch Panel.
H	Height
W	Width
L	Length
Isthmus AP	Used to measure the Isthmus height distance.
Cortical Thickness	Cortical thickness of the lymph node.
Delete	Press to delete this anatomy.

TI-RADS (ACR)

The ACR (American College of Radiology) TI-RADS (Thyroid Imaging, Reporting and Data System) is designed to help the physician in the assessment and systematic documentation of Thyroid nodules based on the peer reviewed published literature, ACR Thyroid Imaging, Reporting and Data System (TI-RADS): White Paper of the ACR TI-RADS Committee by Franklin N. Tessler, MD, CM, Jenny K. Hoang, MBBS, William D. Middleton, MD, Edward G. Grant, MD, Michael D. Beland, MD, Lincoln L. Berland, MD, Sharlene A. Teefey, MD, John J. Cronan, MD, Ulrike M. Hamper, MD, Terry S. Desser, MD, Mary C. Frates, MD, Lynwood

W. Hammers, DO, A. Thomas Stavros, MD, Jill E. Langer, MD, Carl C. Reading, MD, Leslie M. Scoutt, MD, A. Thomas Stavros, MD published on Journal of the American College of Radiology Volume 14, Issue 5, May 2017, Pages 587-595.

The ACR TI-RADS uses standard terminology to document the features of the nodules. Points are given for all the ultrasound features in a nodule, with more suspicious features being awarded additional points. When assessing a nodule, the user selects one feature from each of the first four categories and all the features that apply from the final category and sums the points. The point total determines the nodule's ACR TI-RADS level, which ranges from TR1 (benign) to TR5 (high suspicion). No recommended action is provided by the system to the user. In the ACR TI-RADS publication, recommendations for FNA or ultrasound follow-up are based on a nodule's ACR TI-RADS level and its maximum diameter.

TI-RADS (ACR) is Thyroid Imaging Reporting and Data System (from ACR). ACR TI-RADS could provide guidance regarding management of thyroid nodules on the basis of their ultrasound appearance.


NOTE

TI-RADS is a registered trademark of ACR and all rights reserved by ACR.

To activate TI-RADS on the system:

1. The operator selects the appropriate probe and Thyroid preset.
2. The operator selects Nodule in Measure Menu.

Figure 6-68 Select Nodule



3. The operator measures Nodule L*H*W* manually or Auto Contour(HxL)/(HxW).


4. The operator selects TI-RADS(ACR) on Touch Panel. Then TI-RADS(ACR) menu appears. The operator selects each category for the Nodule. All points will be summed up automatically.

Figure 6-69 TI-RADS(ACR) Menu



5. The operator clicks  icon at the right corner of TI-RADS(ACR) menu, a user prompt window will appear. The window has two display modes.
- Place the cursor on the icon, the window will display a summarized introduction for ACR TI-RADS.

Figure 6-70 Summarized User Prompt Window



Advanced Features

- Move the cursor to the category, the window will display a detailed introduction for the selected category.

Figure 6-71 Detailed User Prompt Window

Category	Points
Cystic or almost completely cystic	0 point
Spongiform	0 point
Mixed cystic and solid	1 point
Solid or almost completely solid	2 points

Total Points

0 Points	2 Points	3 Points	4 to 6 Points	7 Points or More
TR1 Benign	TR2 Not Suspicious	TR3 Mildly Suspicious	TR4 Moderately Suspicious	TR5 Highly Suspicious
No FNA	No FNA	FNA if >= 2.5 cm Follow if >= 1.5 cm	FNA if >= 1.5 cm Follow if >= 1.0 cm	FNA if >= 1 cm Follow if >= 0.5 cm *Rule to discussion of papillary microcarcinomas for < 5 mm TR5 nodules

Level

TR1 TR2 TR3 TR4 TR5

Comments

The operator clicks  icon again to exit the user prompt window.

6. The operator marks TI-RADS Level TRx manually.

NOTE

The determine definitive diagnosis depends on FNA (fine needle aspiration).

7. After the completion, the operator presses **Report** key to check all the measurements and TI-RADS in the Worksheet.

Worksheet and Summary Worksheets

Worksheets and Summary Worksheets are provided for all documented Thyroid anatomies.

To move to the next page, rotate the **Page Change** rotary beneath the Touch Panel.

NOTE

Only defined features are displayed on the Summary Report. To display the undefined features, select **Show Undefined Features** at the bottom of the Summary Worksheet.

NOTE

To exit back to the previous measurement screen, press **Set** on **Exit** button.

NOTE

To exit back to the scan screen, press **Worksheet/Summary** on the Touch Panel.

Measure Assistant Thyroid (Auto Contour)

You can request that the system trace/outline the border of a nodule using Measure Assistant Thyroid (Auto Contour). You do this by setting the Region of Interest (ROI) around the nodule; the system can then measure the nodule by drawing the contour around it.

To automatically detect the nodule on the display,

1. Press Measure.
2. Press Auto Contour (HxL) or Auto Contour (HxW) .
3. Place the Cursor in the center of the nodule and press Set. Size the ROI around the nodule. Use the Trackball to resize the ROI.
 - To increase the size of the circle, move the Trackball down or right.
 - To decrease the size of the circle, move the Trackball up or left.

NOTE

Include the entire nodule, even if additional surrounding tissue is included.

4. Press Set. A trace appears around the nodule.
5. Size the trace via the Trackball.
6. Press Set. The contour around the nodule is generated.

NOTE

Multiple nodule traces may be generated by the system.

7. Inspect the generated contour for accuracy. If edits are necessary, execute steps 8-10 to edit the contour prior to accepting the measurement. Otherwise, skip to step 11.
8. To cycle through the generated contours, press Cursor Select.
9. To edit the selected contour, move the Trackball to appropriately size the edit region and then press Set.
10. The blue portion of the contour can be edited by moving the Trackball to the portion of the contour you want to edit.

NOTE

The Caliper closest to the cursor enables editing.

11. After you have completed your edits, press P1 to accept the measurement.

Scan Assistant

Scan Assistant Manager

Scan Assistant provides an automated exam script that moves you through an exam step-by-step. This allows you to focus on performing the exam rather than on controlling the system and can help you to increase consistency while reducing keystrokes. The system automatically invokes the correct mode and imaging parameters, advances to the next step in an exam, annotates the image, initiates measurements, and assigns the measurements to the worksheet/report.

Availability

The following additional imaging parameters and preferences are available for use in a Scan Assistant program: CW Doppler, Dual on Freeze, Depth, Color Scale, PW Doppler Scale, PW Sample Volume size, and Flow Model Selection.

You can initiate one or more manual Doppler measurements/calculations.

Body Patterns are available for use during a Scan Assistant program. You can turn a Body Pattern on/off, select a particular Body Pattern graphic, and specify the position of the probe mark on the Body Pattern graphic.

The footswitch can be used with Scan Assistant. You can map Pause/Resume, Previous Step, and Next Step to the footswitch.

The “Always Use Doppler Cursor” preset, available on the Utility --> System --> General page, allows all PW Doppler steps to start with full screen 2D image plus mode cursor. You can specify the Store Order in Scan Assistant to set the Reading Order for the radiologist. The Learn Probe attribute can be set to learn and change the probe for the user in the middle of the exam.

Scan Assistant Definitions

Scan Assistant definitions:

- **Scan Assistant Manager.** Available via the Utility -> Scan Assistant page to import/export Programs created via the Scan Assistant Creator and to assign Programs to a user/exam category.
- **Import.** Used to load Programs created via the Scan Assistant Creator on to the ultrasound system.
- **Export.** Used to move Programs from one ultrasound system to another ultrasound system.
- **Scan Assistant Creator.** Used to create Scan Assistant Programs.

Scan Assistant Description

Figure 6-72 Scan Assistant Display Description



1. Program name.
2. Completed steps/out of total number of steps, and step description area.
3. Program step status (Complete/Incomplete), step number, step name. A checkmark indicates that this step has been completed. You can also manually check the box to bypass this step.
4. This column indicates that the action moves the Program to the next step.
5. Active step.
6. This column indicates the mode or when a measurement needs to be made.
7. Edit (Pencil Icon).
8. Stop. Stop also allows the program to be stopped, restarted, or a new program selected.
9. Pause/Resume.

Setting up Scan Assistant


To set up Scan Assistant,

1. Import the Scan Assistant Program created using the Scan Assistant Creator or exported from another ultrasound system program.
 - a. Insert the media with the saved Program from the Scan Assistant Creator or exported program from another ultrasound system.
 - b. Press **Utility > Scan Assistant**.
 - c. Select Import from the **Scan Assistant Manager** page.

Advanced Features

- d. In the Source field at the top of the Import Programs pop-up, select the media that the Program is stored on.


Figure 6-73 Import Programs



- e. Select the Program(s) to be imported. If a folder is selected, all programs in the folder will be imported.

- f. Select Import. The Program(s) you selected are stored to the ultrasound system. Select **Exit**. Then you can add it to the exam category and user.

Figure 6-74 Import Completed




- Assign the imported Program to the exam category and user. Under Program Selections on the right-hand side of the Scan Assistant Manager page, specify the Exam Category and User for this Program. You can select All Users, or a specific user. If you specify All Users, all users will have the ability to use this Program while in the specified exam category, unless the user has his/her own list defined.

Advanced Features

3. Select the imported Program from **Available Programs > Custom Programs** on the left-hand side of the page. Then press the right arrow button to move the imported Program to the exam category and user selected above.

Figure 6-75 Add Program



You can access the Scan Assistant Creator to edit the exam's program from the imaging display via the Creator Icon located at the bottom of the Scan Assistant Program monitor on the display. You can activate the Scan Assistant Creator from the image screen, make edits, and then run Scan Assistant to test your changes.

NOTE

If you edit the program after you have already stored several images, and your edits change the number of program steps, you are prompted to Restart or Continue the Scan Assistant program.

NOTE


If you edit the program after you have already completed several steps, the check marked steps will remain checked, even if you insert a new step in between the currently checked steps. If this is not correct, you can edit the check marks or restart the program.

Using Scan Assistant

After you have set up Scan Assistant, the Program is active when you exit the Patient menu. The Program is located on the left-hand side of the display and as you can see in the

example below, the annotation for the first step has been automatically noted on the image, ready for you to scan the specified anatomy.

Figure 6-76 Scan Assistant Display



1. Follow the steps indicated in the Program: image/measure the appropriate anatomy.
2. Perform the indicated trigger to move to the next step in the Program.

NOTE


The footswitch can be used with Scan Assistant. You can map Pause/Resume, Previous Step, and Next Step to the footswitch.

3. To pause or unpause Scan Assistant, press the pause button on the display.
4. To stop or restart a Program, press the Stop icon at the bottom of the Scan Assistant Program. A dialog pops up. This dialog lets you restart the current Program, start another Program, or stop Scan Assistant.
5. To skip a step or move to a certain step, press the up/down arrows on the keyboard or select the step you want to move to using the Trackball and Set keys.

Activating the Scan Assistant Creator

To activate the Scan Assistant Creator on the ultrasound system, press the **Creator** key.

Figure 6-77 Scan Assistant Manager



Exporting Scan Assistant Programs


Exporting Scan Assistant Programs allows them to be imported to another Versana Premier/ Versana Premier Lotus or to be edited offline with the Scan Assistant Creator tool.

Follow below steps to export a Program.

1. Insert the media to save the program.
2. Press **Utility > Scan Assistant**.
3. Select **Export** from the **Scan Assistant Manager** page.
4. In the Source field at the top of the Export Programs pop-up, select the media that the Program is to be stored on.
5. Specify the Program Directory using the drop-down menu if the desired Program Directory already exists on the media. If not, or if you want to export the Program to a new Program Directory, type a new Program Directory name in the field.
6. Select the Program(s) to be exported. If a folder is selected, all programs in the folder will be exported.

7. Select **Export**. The Program(s) you selected are stored to the media. You can now import it to a new ultrasound system.

Figure 6-78 Export Program(s)



My Trainer

My Trainer provides a quick guide to operate the system.


Follow below steps to access My Trainer:

- Press **Alt + H** to enter **My Trainer**.

OR


- Press **Assistant** key on the control panel, then select **My Trainer**.

Figure 6-79 Enter My Trainer by Assistant key



There are five sections in My Trainer. The five sections are displayed on the left side of the My Trainer interface.

Figure 6-80 Sections Display




1. Getting Started: Supply Mains Switch, Probe Connection.
2. System Setting: Facility name set up, Date time set up, Enable options check, Language setting, User configuration, User preset.
3. Peripheral Connection: Back panel ports, Wire Connectivity, Wireless Connectivity, Connect to a printer, DICOM worklist setting, DICOM Image storage, Back-up/restore the animal information, software version.
4. Maintenance: Trackball, Contact GE HealthCare service, Connect to Insite, Log export, Air Filter Replacement.
5. New Feature: Latest new features.

My Trainer User Interface

The interface of the My Trainer includes below information:

Figure 6-81 Interface Illustration



1. The product name of the system.
2. Sections. Select the button to expand the section to show the list of subsections.
3. Subsections.
4. The highlighted subsection shows this subsection is currently opened and displayed on the right side.
5. Exit. Press this button to exit this interface
6. The left number shows the current page of the subsection. The right number shows the total pages of the subsection.
7. Illustration with graphic.
8. Text instruction about how to operate.
9. Step Description.

Follow-up Tool

The Follow-up tool is intended to perform serial scans on a patient, and compare the images of a previous ultrasound exam to the current exam.

Follow-up tool can select up to 4 historical data for comparison.

When performing a follow-up exam, the system automatically reloads the scanning parameters including ROI for QAnalysis from previous exam, and allows for side by side scanning for image comparison. This allows the physician to use consistent scanning parameters and Region of Interest (ROI) from exam to exam on the same patient and may assist in assessing a patient condition over time.


Follow-up tool configuration

- Follow-up tool is available in B / B (CHI) / CF/ PDI mode.
- In order to utilize the follow-up tool, images should be saved in Raw DICOM format.

To configure image format as Raw DICOM:


Enter **Utility > Connectivity > Button** to set the Print button for this format. Raw Dicom should be selected in the **Format** field.

Figure 6-82 Set Image format




- The Customer can configure the Follow Up Button in **Utility > System > User Configurable Key**.

Figure 6-83 Follow Up Tool Configuration



- The Customer can configure the Follow Up Tool settings in **Utility > System > System Imaging > Follow Up Tool**.

Figure 6-84 Follow Up Tool Setting



- Comparison Image Side: select the comparison image displayed position. The Choices are: Right and Left. The Factory Default value is Left.
- Comparison Date: Choose the date display condition of images which are used for comparison. The choices are: All Dates, Different Dates, None.
 All Dates - The date is always displayed on the comparison image.
 Different Dates - The date only displays when the date of the comparison image is different from the active exam date.
 None - No Comparison image date is displayed.
- Show label on Clipboard: choose when to show the comment label for the image on clipboard. The Choices are: Always, Follow up tool and Off.
- Auto Copy QA Measurement:
 Check the box, QAnalysis measurement is copied automatically from comparison image to the current scanning image.
 Uncheck the box, QAnalysis measurement is not copied automatically from comparison image to the current scanning image.

The factory default is that the box is unchecked.

Follow-up Tool Workflow

The First Exam

1. Press **Patient** to enter the Patient page.
2. Create a New Patient, and create a new exam for the patient. Select **Scan** to enter scanning mode.
3. Press **Probe** to select the application and preset.
4. Scanning with B mode and adjust the scanning parameters.
5. Press **Freeze**.

NOTE

Comments and body mark can be added, and will be stored with the images.

6. Activate CF or PDI mode and adjust the scanning parameters. Here CF mode is used for example.
7. Press **Freeze**. Press the Print key to save the image.

NOTE

Comments and body mark can be added, and will be stored with the images.

8. Press **Measure** to activate measurement and select measurement from measurement folder.
9. Select measurement folder and select a measurement item.

NOTE

Follow-up Tool only available for the measurement items under MSK -> Knee -> Fusion.

10. Draw only one ROI and the result is displayed in result window.
Select **Rt** (Right) or **Lt** (Left) for the patient's right and left side.
Select **Prox** (Proximal), **Mid** (Middle) or **Dist** (Distal) if the vessel has a location.
11. Press the Print key which is assigned to save the image with measurement.

The Follow-up Exam

NOTE

Before starting Follow-up exam, define a key from User Defined Key for activating/deactivating Fusion via **Utility > System > User Configurable Key**.

NOTE

Before starting Follow-up exam, set the Fusion values via **Utility > Imaging > B/HAR** mode page.

1. Enter Patient page, and select the patient from the **Patient View** list.
2. There are different methods to activate Follow up exam.

NOTE

Be sure to activate the same probe used for the initial exam.

- To access Follow up exam:
 - a Select the Exam.
 - b Select **Image History** or **Active Images**.
 - c Select one image, then select **Compare**.

NOTE

Only one image is selected at a time so **Compare** can be activated. If more than one image is selected, **Compare** can not be activated.

- In scanning mode, select the Follow up exam icon, the Follow up exam is activated.

Figure 6-85 Follow up exam icon



OR

Press **Follow Up** key on the control panel to activate Follow up exam.

3. The system enter scanning mode with Follow-up exam. The screen reverts to split screen.

Click the Follow-up icon to deactivate or activate the Follow-up exam.

As factory default setting, the live image is displayed on the right side of the screen and the selected comparison image is on the left.

The user can configure the display position in **Utility > System Imaging > Follow Up Tool > Comparison Image Side**.

4. Select a B mode image on the clipboard to do follow up and the body mark and the comments which are added in the selected image will be transferred to the current scanning image.
5. Press the key assigned for Fusion to turn on Fusion. The overlay will appear on the active image. Align the bony landmarks or fascial planes with overlay to find the same scanning position from comparative study.

NOTE

You can assigned a key from key User Defined Key to turn on/off Fusion via **Utility > System > User Configurable Key**.

Move the probe until the scanning position gets to the mark.

6. Select the assigned key again to deactivate Fusion.

NOTE

Fusion should be turned off before saving the image, otherwise the mark will remain on the image which is frozen and saved.

NOTE

For additional follow-up, select the appropriate clipboard image for comparison.

7. Press **Freeze**. Press the Print key to save image. An image with split screens will be stored in the exam report.
8. Select one image of Color Flow Mode on the clipboard to do follow up. Color Flow Mode is used for example here.

The color window automatically appears on the same position and in the same size with the comparison image.

The scanning parameters is reloaded from the comparison image to the current image.

9. Select Fusion from the assigned Function key. When Fusion is activated with image that has been saved in Color, the color disappears automatically from color window. The color bar and color window will remain.
10. After the scanning position is located, press the key assigned for Fusion to turn off Fusion. When the Fusion is turned off, color appears automatically in color window. During this operation, the user does not need to press **Color** key to activate color.
11. Press **Freeze**.
Select the desired frame by using **Frame By Frame** control on the touch panel.
12. Press **Measure**. The measurement method, measurement item and ROI will be transferred automatically from comparison image to the current scanning image.

NOTE

Be sure Auto Copy QA Measurement (Dual Screen) has been checked in the Follow Up Tool portion in **Utility > System > System Imaging**. This is not available through factory default.

NOTE

If Auto Copy QA Measurement (Dual Screen) is not selected, the QAnalysis measurement will not be copied automatically from comparison image to the current scanning image. The user can perform the measurement with new measurement method, measurement item and ROI on the current image.

The result window displays the measurement result of comparison image and the current image.

13. Press the Print key to save the image and measurement.

The image with split screens and measurement result for current image and comparison image is stored.

Worksheet

To view a worksheet, press the **Report** key on the Touch Panel.

Report

To view the report, press **Report** on the Touch Panel.

Select **Whizz Report** on the Touch Panel to enter report. Select the image and click the + icon, the image will display in the **Images** section.

The measurement result window with measurement data of current image and comparison image is kept on the image and displays.

QAnalysis (Option)

Quantitative Analysis (QAnalysis) is available for the image loop acquired in the following modes: Color Flow Mode, Power Doppler Mode, Tissue Velocity Mode, Strain and Strain Rate. The difference of QAnalysis operations is small among different modes.

Activating QAnalysis

1. Scan the patient in CF/PDI/TVI mode or select a desired cine loop in CF/PDI/TVI mode from the stored images.

NOTE

Images from the current scan session acquired (already in CINE) or from a saved image loop can be used for QAnalysis.

NOTE

QAnalysis is only available if the user has selected an image loop. If the user has selected a saved still image (just one frame), QAnalysis is not available.

2. Press **Freeze**. Then move the trackball to enter CINE mode.

NOTE

QAnalysis is only available when the system is in CINE mode.

3. Select **QAnalysis** to activate Quantitative Analysis function. To toggle the trackball function between QA and Scroll, press the **Set** key.

QAnalysis Screen Description

Figure 6-86 TVI mode QAnalysis Example




Table 6-27 QAnalysis Screen Description

1.	TVI Cineloop Window Sample Area: Indicates sampling position of the velocity. The sample area is color-coded: the first sample area is yellow, the second green, etc.
2.	B Cineloop Window Sample Area: Indicates sampling position of the velocity. The sample area is color-coded: the first sample area is yellow, the second green, etc.
3.	Analysis Window. - Y axis: Velocity scale (cm/s) - X axis: Time(s) - Time at cursor position. - Velocity at Cursor position. - Velocity at frame marker position (Color coded)
4.	Sample Area
5.	Time at cursor position and velocity at cursor position. Position the pointer cursor over the analysis window.

6.	Sample Area Tools. - Pencil Icon: Creates a sample area based on freehand drawing. - Shape Icon: Creates a sample area with a pre-defined circular/ellipse shape.
----	---

Figure 6-87 CF Mode QAnalysis Example



Selecting QAnalysis Image Range

A range of frames is selected for the QAnalysis in Cine mode (before accessing QAnalysis). Only the frames in this range are used for the QAnalysis.

If a range is not selected prior to accessing the QAnalysis, the system uses the default Cine start and end frames as the default start and stop frames.

1. The first frame in the analysis series is selected by adjusting the **Start Frame** control to the desired frame
OR
using the **Trackball** or the **Frame by Frame** control to select the desired first frame and then selecting the **Start Frame** control.
2. The last frame in the analysis series is selected by adjusting the CINE **Last Frame** control to the desired frame
OR
using the **Trackball** or the **Frame by Frame** control to select the desired last frame and then selecting the **Last Frame** control.

Generating a Trace

Trace from a pre-defined sample area

1. Select the sample area Ellipse ROI button (shape icon on the monitor display).
2. Move the cursor to one of the Cineloop windows using the **Trackball**.
3. Press **Set** to anchor the sample area.
In this frame, the sample area is marked with an anchor. If the cineloop has more than one heart cycle, a sample area will also be anchored in the corresponding frame in the next heart cycle.

The trace is updated accordingly in the Analysis window.

Trace from freehand sample area

1. Select the Freehand ROI button (pencil icon on the monitor display).
2. Move the cursor to one of the Cineloop windows using the **Trackball**.
3. Trace the outline of the desired ROI by moving the caliper with the **Trackball**.
4. Press **Set** to anchor the sample area.
The sample area is automatically closed and the trace is updated accordingly in the Analysis window.

Manual tracking of the sample area (dynamic anchored sample area)

1. Place a sample area over a region of interest. Note the anatomical location of the sample area.
2. Scroll to a new frame using the **Trackball**.
3. Press **Scan Area** key until the QA trackball assignment is selected.
4. Move the cursor to the sample area using the **Trackball**.
5. Press **Set**. The sample area is unanchored.
6. Drag the sample area to the corresponding anatomical location in the new frame.

When the sample area is anchored in more than one frame, linear interpolation is performed so that the sample area is smoothly moved between the anchored positions in the selected frames when running the cineloop.

7. Press **Scan Area** key until the scroll trackball assignment is selected.
8. Using the **Trackball**, scroll through the cineloop and control that the sample area follows the moving anatomical structure.
9. Add anchored sample areas in several frames to obtain a more accurate displacement of the sample area.

Moving a dynamic anchored sample area

1. Press **Scan Area** key until the scroll trackball assignment is selected.

2. Using the **Trackball**, browse through the cineloop to display one of the frames where the sample area was anchored.

NOTE

In these frames, the sample area is marked with an anchor.

3. Press the **Scan Area** key until the QA trackball assignment is selected.
4. Move the cursor to the sample area using the **Trackball**.
5. Press **Set**. The sample area is unanchored.
6. Drag the sample area to a new location.
7. Press **Set** to anchor the sample area to the new location.

Delete a trace

The user can delete all traces at once or one at a time.

1. Move the cursor over one of the sample area. Confirm that cursor is changed to hand icon.
2. Select **Delete Sample Area**.
3. Select **Current Sample** or **Delete all** as necessary.

Manipulating the Sample Area

Up to eight ROIs can be saved on the reference image, with the corresponding eight traces plotted simultaneously on the graph. Each ROI display has a different color, and its corresponding trace data is plotted using that same color.

Once eight ROIs have been saved, the system does not automatically generate an active ROI when the cursor is positioned over the displayed reference image.

The saved ROIs can be a mixture of elliptical and freehand ROIs.

When the user repositions an ROI, the old trace data is erased from the plot and the trace data for the new position replotted.


If the ROI position on the last frame of the selected image range is moved, the corresponding ROIs on all frames are repositioned to match the last frame.

The user shall also have the capability of setting separate ROI positions on different frames of the contrast images, and the system shall linearly interpolate the ROI positions for the frames in between the selected frames.

Setting the default sample area shape


1. Position the cursor on the ROI on Cineloop windows, a white circle appears.

Figure 6-88 Sample Area Information Box



2. Select **Set sample area shape** on the System menu. The Information Box displays.

Figure 6-89 Sample Area Information Box



3. Select Height, Width and Tilt angle.
4. Select **Set as default**. The current ROI size is set as the default for subsequent Ellipse ROIs.

Reshaping a Sample Area

To reshape the sample area:

1. Position the cursor on the ROI to reshape and press **Set** key.
2. The ROI system menu displays. Select **Set sample area shape**.
3. Adjust Height, Width and Tilt angle.

4. Press **OK**. The selected ROI size changes.

Sample Area Shapes

There are two different methods for determining the shapes of the sample area.

Ellipse ROI

1. Select the ellipse icon (shape icon on the monitor display).
2. When the trackball positions the image display cursor over the reference image(s), an elliptical ROI is automatically generated and displays on the reference image(s).
3. The average velocity value inside the ellipse is calculated for every image in the image analysis range and plotted in the image display area.
4. The last generated or selected ellipse is considered the active ROI, and its trace plot automatically updates as the user repositions it on the reference image. Old traces are erased.
5. When scanning with an elliptical ROI, press **Set** to fix the ROI position and freeze its corresponding trace on the plot. A new active ROI is generated whose position is manipulated by the trackball and whose velocity curve traces will be plotted as before, while the previous ROI and trace remain fixed at the points they were saved at.

NOTE

Elliptical ROIs can be positioned in any manner that keeps their center within the image boundaries. In the case that part of the ROI is outside the image boundary, only data from within the image boundary is used for calculating the mean velocity value.

NOTE

You can change the size of the Ellipse ROI by adjusting the Ellipse control.

Freehand ROI

1. Select Freehand icon (pencil icon on the monitor display).
Use the **Trackball** to position the caliper on the reference image at the start point. Press **Set** to fix the start point.
2. Trace the outline of the desired ROI by moving the caliper with the **Trackball**.
3. When a suitable ROI has been drawn, release the **Set** key.
The caliper is then free for repositioning for another freehand ROI.

NOTE

You cannot go outside the image boundary when drawing a freehand ROI.

Deleting a Sample Area

Sample ROIs and their corresponding traces can be deleted using **Delete Sample Area**.

1. Select **System Menu**; a pull-down menu displays.
2. Select **Delete Sample Area** to delete the currently active ROI.
Select **Delete all Sample Area** to delete all currently set ROIs and all of their traces.

NOTE

The corresponding traces for the deleted ROIs are erased from the plot.

NOTE

Deleting an ROI causes the ROIs to be deleted from all frames in the analysis loop.

Disabling/Enabling the frame

NOTE

Disabling/Enabling the frame is available when enter QAnalysis from CF and PDI mode.

Frame disabling excludes the actual frame from the cineloop display. Frame disabling is available only with contrast data.

Disabling the frame from the frame marker

To disable One Frame:

1. Use the **Trackball** to move the cursor to the frame maker to disable.
2. Press **Set** to disable the frame.
3. The frame marker is changed from green to red to indicate the frame has been disabled.

NOTE

The disabled frame is no longer displayed in the reference window when scrolling through CINE memory.

Disabling multi-frames from the frame marker

1. Use the trackball to move the cursor to the first frame maker to disable.
2. Press and hold down **Set**.
3. Move the cursor with the Trackball to the last frame to be disabled and release Set.

The marker is turn red and the data from that frame is removed from the trace and any subsequent trace processing.

Disabling a frame from the cineloop window

1. Use the trackball to move the cursor to the cineloop window.
2. Select **Disable frame**.
The current frame is disabled and the corresponding frame marker displays red.

To enable the frames

To re-enable all deleted frames:

1. Position the cursor on the Frame Marker line.
2. Select **Enable all frames**.
3. All disabled frames are re-enabled.

Smoothing

NOTE

Smoothing is not available when entering QAnalysis from CF/PDI mode.

The system can smooth the traces displayed by applying a filter over a defined time window. The type of filter available is depending on the analysis signal displayed.

1. Select **Smoothing**.

NOTE

When smoothing is turned on, it applies to all traces in the plot window.

2. The smoothing filter list displays. Select the appropriate parameter.

Horizontal Sweep

Horizontal Sweep allows you to increase or decrease the time interval over which to plot the analysis curve.

The Horizontal Sweep control can range from TBD on the short side to the time interval between the user selected first and last frame. The default is the user selected image range. If the user has not yet selected a first and last frame, the first and last default frames from the displayed CINE loop are used.

Drift Compensation

NOTE

Drift Compensation is only available when entering QAnalysis from TVI mode.

Drift Compensation compensates drifting of Tissue Tracking curves by either resetting the curve to zero at the tracking start point (cycle resetting) or by linear compensation throughout the cycle (linear compensation).

NOTE

When Displacement is chosen by Analysis Signal, Drift Compensation is active.

Statistics

Select **Statistics** to enable/disable display of statistics of the frame or loop. The statistics are shown only when the loop is stopped.


- Ratio: Ratio of Color (Power) Doppler pixel over total ROI area.
- Area (mm²): The size of ROI.
- Max Ratio/ Time of Max Ratio: Maximum Ratio of Color (Power) Doppler pixels in each ROI, and which frame that occurs in.
- Min Ratio/Time of Min Ratio: Minimum Ratio of Color (Power) Doppler pixels in each ROI, and which frame that occurs in.

Trace Measurements

Max Gradient

Displays the time and gradient that becomes the maximum gradient between the CINE start and end frame.

Figure 6-90 Max Gradient




Exporting Traces (Saving the Trace Data)

You can save the trace data to an external file.

1. Select **Export Traces** on the touch panel to save the trace data.

- The following window displays.

Figure 6-91 Export Trace Window



- Location: Select Location which to save.
 - Filename: Enter the filename. (Only Text)
- Select **OK** to save the data and return to the QAnalysis screen.
 - All displayed ROI traces are saved in the exported file.

All plot data (intensity, gradient and gradient derivative) are exported to a text file by “Export Trace”.

Table 6-28 Example of exported file

Time(s):	Trace 1:	Trace 1 dGrad.:	Trace 1 dGrad.
0.00000	-3.97995e+000	-2.15924e+001	8.05159e+001
0.03121	-5.14631e+000	-1.64719e+001	1.74256e+001
0.06242	-5.75798e+000	-1.27675e+001	-7.78004e+001
0.09362	-6.02222e+000	-1.27675e+001	-1.93426e+002
0.12483	-6.11224e+000	-1.44515e+001	-4.17252e+002

NOTE

The Smoothed trace is the one saved if the user has applied a smoothing filter.

NOTE

Only data from the user selected image range is included in the exported trace file.

NOTE

No trace results are saved in the standard image database.

Annotating the QAnalysis Data

The user can annotate both the reference image and the trace plot displays. Use **Comment** key to type the annotation. See Chapter 6 for reference.

Exiting QAnalysis

There are several methods to exit QAnalysis.

- Select **Exit QAnalysis**.
- Press **Freeze** to unfreeze and resume scanning.
- Press any other button that returns the system to real-time scanning.

Whizz Label (Not applicable for China)

The right upper quadrant (RUQ) is an anatomical area where a patient presents with sign or symptoms suggesting clinicians to take look or investigate. Ultrasound is a non-invasive imaging method frequently used to image this anatomical region.

Whizz Label is a software technology utilizing Artificial Intelligence (AI) during ultrasound acquisition to automatically detect the liver, gallbladder, and/or right kidney during an ultrasound exam.

This technology is useful for clinicians that need to label or annotate ultrasound images for documentation purposes similar to manual annotation of ultrasound images today.

NOTE

Whizz Label should not be used for pediatric abdomen scanning. It is only supported for Adult abdomen scans.

NOTE

Algorithm cannot detect proper anatomy when probe movement is too fast or unstable. Algorithm will not identify labels if the image quality is not deemed acceptable, e.g. shadows, excessive or very little gain. In addition, presence of pathologies or missing parts of the kidney capsule may impact auto label performance.

NOTE

Whizz Label is only available for the RUQ organs. Accuracy of figures is specified in the context of right upper quadrant images and it is not applicable to the scanning of other parts of the body.

NOTE

To use Whizz Label, the right upper quadrant organs are expected to be imaged as follows:

- Liver - in epigastric right midclavicular line or right intercostal windows using longitudinal (sagittal), transverse views and/or oblique views.
- Gallbladder - scanned from subcostal or right intercostal approach.
- Right Kidney - scanned in epigastric/subcostal or right intercostal windows using longitudinal (sagittal), transverse, and/or oblique views.

The tool has not been tested when scan is performed from other scan windows. Specifically, the algorithm will not perform as expected when scanning the abdomen RUQ from a subxiphoid window

NOTE

Whizz Label needs to see large portion of the organ to be accurate, as per the size requirements shown in the table below.

Table 6-29 Accuracy of Whizz Label

Organ	Accuracy
Liver	90% Accuracy when > 32cm ² size section of the liver is visualized on the image scan area.
Gallbladder	80% Accuracy when > 8cm ² of the longitudinal view of gallbladder detected.
Right Kidney	90% Accuracy when > 16cm ² size section of the right kidney is visualized inferior to liver in longitudinal or transverse view(s).

Using Whizz Label

User can start Whizz Label to label liver, gallbladder and/or right kidney in frozen image.


Whizz Label supports multi-language but can only be used in the Abdomen preset with specific probes in B/CHI/CF/PDI mode. Refer to *Table 5-8 Probe modes of operation and Features* on page 330 for detail information.

Whizz Label (Freeze)

To use Whizz Label on a Frozen Frame:

1. Select Whizz Label on Touch Panel. Once image is frozen, system will label liver, gallbladder or right kidney as seen on the frozen image automatically.


Figure 6-92 Whizz Label while Freezing



Advanced Features

2. If you are satisfied with the labels, press P1 to save. A confirmation window may display if you configured in **Utility > System > System Imaging > Need confirmation before storing image under Whizz Label**.
You can also activate Comment to modify or move the labels as needed.

Figure 6-93 Whizz Label confirmation



Whizz Label (Live)


To use Whizz Label in Scan Coach while scanning, user needs to enable the function in **Utility > System > System Imaging > Enable Whizz Label in Scan Coach**.

NOTE

When storing a cine loop in Scan Coach, auto generated labels & protocol label will not be stored.

1. Once Whizz Label in Live mode is set to ON in the configuration, organ labels are automatically displayed during scanning when in the relevant protocol step (for Liver, Right Kidney or Gallbladder). When the user Freezes a scan, the labels will be shown only for the last frame.

Figure 6-94 Whizz Label while Scanning




2. If user is satisfied with the labels, user can press P1 to store the frozen image. The system will show a message window to ask user to confirm correctness of labeling

before image storage. Whizz Label will be turned off when saving the image. System will resume Whizz Label after storing

User can configure whether to display the confirmation window in **Utility > System > System Imaging > Need confirmation before storing image under Whizz Label**.

You can also activate Comment to modify the labels as needed.

Figure 6-95 Whizz Label confirmation



3. If you switch modes during scanning or turn on crossbeam, system will turn off Whizz Label and remove the labels.

If you enter comment or body pattern during scanning, system will turn off Whizz Label but keep the labels, user can also modify the labels kept on the image.

Advanced urology procedures

E7C8L-RS probe can be attached to a stepper to be used for urological procedures including transperineal needle guidance, transperineal grid biopsy and prostate brachytherapy.

ERB7/ERB Biopsy Attachment for E7C8L-RS


ERB7/ERB Biopsy attachment can be installed on the E7C8L-RS probe to assist you in conducting biopsy procedures. Using this attachment allows a biopsy position to be selected along the guide line while monitoring the Linear and Convex images in real-time.

ERB7/ERB Biopsy Attachment User Instruction

For ERB7/ERB Biopsy Attachment user instruction, please refer to Common Reusable Biopsy Attachment User Manual for more information.

Please pay attention to the alignment tips below:

Figure 6-96 Probe/Attachment Alignment



a Do not install the guide attachment out of the scale of the probe.

b Align the scale of the probe and guide attachment.

NOTE

Attachment position can be adjusted based on patient size or distance to anatomy.

For proper setup recommendation for E7C8L-RS grid biopsy, please visit CIVCO website:
<https://www.civco.com/catalog/> for User Instruction.

NOTE

Make sure probe and grid are attached correctly and tightly on the stepper.

A typical exam consists of:

The needle placement guide should match with the electronic needle placement grid displayed on the Ultrasound system monitor. Grid is optional for stack of coins measurements. When used with the stepper volume calculation option, a “stack-of-coins” method can be used to calculate the volume of the prostate gland.

- Placing the patient in the supine or dorsal lithotomy position.

- Carefully inserting a properly prepared E7C8L-RS probe and orienting it to scan the prostate at the base position.
- Position the mechanical stepping device and secure the E7C8L-RS probe to the stepper.
- Make as many stepper volume measurements as necessary to calculate the volume of the prostate.
- Retract the probe (utilizing the stepper device) the necessary step between each incremental slice of the gland.
- Once the entire gland has been evaluated detach the probe from the mechanical stepper and carefully remove it from the patient.

Biopsy Guideline Type Preset Selection

A preset, provided in Biopsy Kit on touch panel, is available for choosing the Biopsy Guideline type.

NOTE

Changing the needle guide type may be duplicated by programming a User Define key. When the User Define key is pressed, the programmed needle guide type will be changed with single key operation.

Figure 6-97 Setting the Biopsy Guideline Type for E7C8L-RS (-C: Convex Transducer)





Figure 6-98 Setting the Biopsy Guideline Type for E7C8L-RS (-L: Linear Transducer)



E7C8L-RS Grid Biopsy Accessories



WARNING

The grid is single use only. Please choose the specific grid for your clinical usage.

The E7C8L-RS biopsy guide line accuracy absolute error shall be $\leq \pm 1.5\text{mm}$, needle guidance accuracy for grid shall be not larger than 1.8mm when using stabilizer/stepper for procedure.

The table below lists the accessories included in E7C8L-RS MPWs Stepper Grid Biopsy.

Table 6-30 E7C8L-RS Grid Biopsy Accessories - MPWs Stepper

Multi-Purpose Workstation Stepper Grid Biopsy Accessory	CIVCO Part No.	Graphic
Multi-Purpose Workstation Stepper	642-510	
Multi-Purpose Workstation Stepper - Cradle Only kit	642-509	
Sterile Grid (Disposable - one time use only)	610-905 (17GA)	
	610-906 (18GA)	
Multi-Purpose Workstation with adjustable floor stand	610-974	

The table below lists the accessories included in E7C8L-RS Classic Stepper Grid Biopsy.

Table 6-31 E7C8L-RS Grid Biopsy Accessories - Classic Stepper

Classic Stepper Grid Biopsy Accessory	CIVCO Part No.	Graphic
Classic Stepper	642-512	
Classic Cradle Only kit	642-511	
Sterile Grid (Disposable - one time use only)	610-905 (17GA)	
	610-906 (18GA)	
Micro-Touch adjustable dual-sided table mount	610-911	
Micro-Touch transportation stand (dual-sided table mount)	610-911S	

The table below lists the accessories included in E7C8L-RS EX3 Stepper Grid Biopsy.

Table 6-32 E7C8L-RS Grid Biopsy Accessories - EX3 Stepper

EX3 Stepper Grid Biopsy Accessory	CIVCO Part No.	Graphic
EX3 Stepper	642-514	

EX3 Stepper Grid Biopsy Accessory	CIVCO Part No.	Graphic
EX3 Cradle Only kit	642-513	
Sterile Grid (Disposable - one time use only)	610-905 (17GA)	
	610-906 (18GA)	
Micro-Touch adjustable dual-sided table mount	610-911	
Micro-Touch transportation stand (dual-sided table mount)	610-911S	

Calibration Procedure

For installing the E7C8L-RS probe on the Grid Biopsy guideline, please visit CIVCO website: <https://www.civco.com/catalog/> for User Instruction.

NOTE

Follow cleaning and disinfecting process for probe, steppers and reusable accessories properly to avoid disease transmission between patients.



WARNING

The biopsy location/position must be calibrated prior to performing the actual biopsy to compensate dimensional errors on kit and probe.



WARNING

Please do grid calibration step by step and check the accuracy before grid biopsy examination.

Depth Calibration

1. Submerge the probe and E7C8L-RS Grid Biopsy in 47°C (116.6 °F) water.

Figure 6-99 E7C8L-RS Grid Biopsy in water




2. Insert the 18 gauge needle into G11 of the grid biopsy.

NOTE

Please make sure choosing the correct needle according to the grid used with grid biopsy kit.

3. Using the right/left knob of the grid biopsy bracket, adjust the location of the needle image corresponding to the location of the grid.


Figure 6-100 Needle Image



Center & Symmetry Calibration


1. Insert two needles on the E7C8L-RS grid biopsy guide in the following locations: B6 and L6.

Figure 6-101 Center Calibration I



1. Relocate the probe on the bracket to adjust the center as shown below. Once it is completed, tighten the knob to fix the location.
2. User can also use Alt + F7 Key and adjust the center by the Left/Right/Up/Down key on the alphanumeric keyboard. The grid displayed on the monitor will appear red indicating it is adjustable. Once position is set, use Alt+F7 again to accept.
3. Repeat the process until the location of the needle exactly coincides with the guide mark.

Figure 6-102 Needle Image



Stepper Volume Measurement

Stepper volume (STVOL) is the method used to calculate the volume of an organ using urology software, a transaxial probe and a mechanical stepping device that moves the probe in fixed increments.

The volumes between the slices are assumed to be segments of a paraboloid of revolution, where the volume between any slice n and $n+1$ is $(A_n + A_{n+1})d/2$.

The sum of the volumes between slices is:

$$V = d \left[\frac{(A_1 + A_2)}{2} + \frac{(A_2 + A_3)}{2} + \dots + \frac{(A_{n-1} + A_n)}{2} \right]$$

Therefore, the total volume of the organ is calculated as:

$$V = d \left[\frac{5}{6}A_1 + A_2 + A_3 + \dots + \frac{5}{6}A_n \right]$$


d : stepper volume increment;

V : stepper volume;

A : area.

Calculation method: Ellipse, Trace, Caliper, Circle, Spline.

Figure 6-103 Stepper Volume Measurement




Stepper volume (STVOL) is the method used to calculate the volume of an organ using the urology software.

To perform a stepper volume (STVOL) measurement:

1. Perform the scan in the appropriate scan plane (Convex plane).
2. Press **Freeze > Measure**.
3. Select the **STVOL** folder on the Touch Panel, an active caliper displays.

4. Perform a STVOL area measurement. For the detail of area measurement, please refer to *Circumference Measurement (Open Trace)* on page 97.
5. The system displays the STVOL area and calculated STVOL in the Results Window.

Figure 6-104 STVOL Area Measurement



6. Press Freeze to unfreeze the image.


Repeat steps 1, 2, 5 and 6. To finish the Stepper Volume measurement, select STVOL Done.



WARNING

Do not freeze the image during a procedure, including biopsy procedure.

Figure 6-105 STVOL Done



NOTE

The maximum number of stepper volume slices is 50.

7. If it is required to perform the STVOL measurement again, repeat steps 1 through 7.
The relative error accuracy for B mode stepper volume measurement is $\leq 30\%$.

To view the stepper volume worksheet, press Worksheet on Touch Panel.

Figure 6-106 STVOL Worksheet


PSA		0.000ng/ml		PPSA Coef(1)		0.12		PPSA Coef(2)		0.12	
Parameter	Value	m1	m2	m3	m4	m5	m6	Method			
B Mode Measurements											
☒ STVOL	2	ml	2.39								Avg.

Strain Imaging

Strain calculates and color-codes the extent of tissue deformation (lengthening or shortening) relative to the original size over a given time interval, typically the systole. Users perform an assessment of myocardial function based on color-coded images.

Strain is calculated using the following standard formula to produce a color map to show the approximate values. The purpose of this tool is not to provide the exact strain or strain rate value.


Figure 6-107 Standard formula for Strain



Advanced Features


Strain is dimensionless index of change in length. It defined as a measure of the amount of local deformation of an object caused by an applied force. It displays the percent of deformation [%].

Figure 6-108 The Strain mode screen



1. Strain color bar
2. Strain start and end markers
3. Strain start and end from detected QRS and Strain sample size
4. Parameter window

Figure 6-109 The Strain Touch panel



Using Strain

Connect and configure ECG before using Strain.

1. From **TVI** mode, press **Strain** on the Touch Panel.
2. Adjust **Strain Start** close to the R-peak.
3. Adjust **Strain End** to end systole, typically near the T-wave.
4. Use the trackball key to adjust the size and position of ROI frame over the area to be examined.

NOTE

User can select a desired cine loop in **Strain** mode from the stored images and enter **QAnalysis** from the Touch Panel

Optimizing Strain

- From an optimized Strain rate display adjust strain tracking to pick out the systolic phase.
- The main use of Strain is to map negative systolic deformation. This means that Strain start and Strain end should be adjusted to pick out the systolic phase of the cardiac cycle: Adjust Strain start close to the R-Peak. Adjust Strain end to end systole, typically near the T-wave.
- The maximum deformation that is color-coded can be adjusted using the Strain scale control. If set too low, most of the wall will show the color indicating maximum deformation. If set too high, the maximum deformation color is never reached.
- Strain provides information only in the beam direction. The apical view typically provides the best window since the beams are then approximately aligned to the longitudinal direction of the myocardium (except near the apex).
- Low strain values may be masked out with a different color using the SI Reject control.

Strain rate

Strain rate calculates and color-codes the deformation per system time i.e. the speed at which the tissue deformation occurs. Users assess myocardial function based on strain rate color-coded images.

Strain rate is calculated as the spatial gradient of velocity data.

Strain rate is dimensionless index of change in length. It defined as a measure of the amount of local deformation of an object caused by an applied force. It displays the speed of deformation [1/sec].


Figure 6-110 The Strain rate mode screen



1. Strain rate color bar
2. Strain length and Strain rate reject

3. Parameter window

Figure 6-111 The Strain rate Touch panel



Using Strain rate

1. From **TVI** Mode, press **Strain Rate**.
2. Use the trackball key to adjust the size and position of ROI frame over the area to be examined.
3. Use the trackball to adjust the dimension of the ROI.

NOTE

User can select a desired cine loop in **Strain Rate** mode from the stored images and enter **QAnalysis** from the Touch Panel

Optimizing Strain rate

- To reduce quantification noise (variance), the Nyquist limit should be as low as possible, without creating aliasing. To reduce the Nyquist limit, reduce the scale while in TVI.
- To check for aliasing, freeze the loop and apply velocity trace (Press Freeze and Q-Analysis), see 'Quantitative analysis'.
- Strain rate provides information only in the beam direction. The apical view typically provides the best window since the beams are then approximately aligned to the longitudinal direction of the myocardium (except near the apex).
- There is a trade-off between noise and spatial resolution controlled by the Strain length control. To minimize noise the Strain length should be maximized.
- The maximum Strain rate that is color-coded can be adjusted using the SRI scale control. If set too low, most of the wall will show the color indicating maximum Strain rate. If set too high, the maximum Strain rate color is never reached.
- Low strain rates may be masked out with a green color using the SRI Reject control.

Using 4D



WARNING

Please take extra care on the volume probe surface temperature before and during scanning to avoid heat injury.

When using volume probe for 4D acquiring, end the 4D acquiring once the current exam is complete, or freeze the 4D image when the machine is not in use.

DO NOT keep volume probe working in 4D acquiring mode for a long time, or the probe surface temperature will exceed the normal value.

4D provides continuous, high volume acquisition of 3D images. 4D adds the dimension of “movement” to a 3D image by providing continuous, real-time displays. With 4D, you can apply rendering techniques to smooth out the appearance of an anatomical structure, for example, a baby’s spine.

You can perform the following types of volume acquisitions within the 4D feature:

Table 6-33 4D Package Options

4D Type	Description	Acquisition Mode
4D	Designed for continuous volume acquisition of a 3D image.	B, 4D
Static 3D	Designed for single volume acquisition of a 3D image.	B, 3D

Features supported with 4D

The following features are supported with 4D:

- Most B-Mode controls
- Annotations
- High Detection Speckle Reduction Imaging (SRI-HD)
- Measurements and Calculations

The following post-processing controls are available with 4D:


- CINE
- Zoom

4D Principles of Operation

The acquisition of volume starts with a 2D image using special probes designed for performing 3D sweeps and 4D scans. The volume box defines the region of interest to be used for the volume sweep.

Volume sweep refers to the range of the sweep of the 2D image to be transformed into a rendered, 3D or 4D image. Static 3D acquisition involves a single volume sweep. 4D involves multiple, continuous volume sweeps.

Figure 6-112 Volume Sweep



- 1 Central 2D Scan
- 2 Start 2D Scan
- 3 Range of VOI Sweep

When initiating a volume sweep, the angle of the volume can be adjusted.

Interactive 3D Rendering


Interactive 3D rendering allows visualization of certain structures, helping to view and analyze different sections of the volume.

Region of Interest (ROI)/Render Box

The Region of Interest (ROI) - also referred to as the Render Box in rendering - contains the section of the volume you want to render. Therefore, objects that are not inside of the box are not included in the render process and are cut out (this is important in surface mode to allow a free line of sight). This may or might not be the entire Volume of Interest (VOI).

The view direction of the ROI can be adjusted.

Figure 6-113 Render View Directions



- | | |
|---|------------|
| 1 | Up/Down |
| 2 | Left/Right |
| 3 | Front/Back |
| 4 | Down/Up |
| 5 | Right/Left |
| 6 | Back/Front |


Render View

In Render view, only the rendered image displays - no reference images.

Image Orientation

Orientation of Image in Sectional view.

Figure 6-114 Quad Render Visualization Mode




- 1 Longitudinal
- 2 Transverse
- 3 Coronal

Principle of Sectional Planes

Sectional planes represent three different planes of the same 3D volume. There are three separate planes, A (Longitudinal), B, (Transverse) and C (Coronal).

Figure 6-115 Illustration of Sectional Planes



The presentation of three orthogonal sectional planes is different from the conventional patient orientation in 2D sonography.


NOTE

Whenever the usual, longitudinal section of the patient to display in field A is selected, the conventional orientation for longitudinal and transverse sections is valid.

Reference Images

Reference images are the individual image displays within the corresponding sectional plane. Reference image A represents the longitudinal view; reference image B the transverse view, and reference image C represents the coronal view.

Figure 6-116 Monitor Display of Reference Images in Sectional View




1. Reference Image A (Longitudinal)
2. Reference Image B (Transverse)
3. Reference Image C (Coronal)

Orientation Help


When you view a 4D image on the display, it's sometimes difficult to recognize the orientation. To help, the system displays a three-dimensional drawing to illustrate the orientation. This drawing displays ONLY in sectional view.

Figure 6-117 Reference Image A




For Reference image A, the transducer plane migrates from the FRONT to the REAR through the volume body.

Figure 6-118 Reference Image B



For Reference image B, the transducer plane migrates from the LEFT to the RIGHT through the volume body.


Figure 6-119 Reference Image C



For Reference image C, the transducer plane migrates from the TOP to the BOTTOM through the volume body.

Examples of Probe Orientation with Reference Planes

Figure 6-120 Abdominal Probe Orientation




4D Operational Controls

Control Panel Overview

When entering 3D/4D mode, the behavior of some of the Control Panel buttons changes. For example, in 3D/4D mode, the PW-Mode, CF-Mode and M-Mode buttons are used (along with Depth) to manipulate the Volume of Interest (VOI).

Figure 6-121 Control Panel Buttons



1. **M** key, M-Mode key used to rotate about the X axis.
2. **PW** key, PW-Mode key used to rotate about the Y axis.
3. **CF** Key, CF-Mode key used to rotate about the Z axis.
4. **3D/4D** key.
5. Trackball, used to move the VOI. Also, the keys surrounding the Trackball map to additional functionality, as shown on the monitor display.
6. **Freeze** key, used to freeze a 4D image.

4D Monitor Display

Imaging parameters are displayed in the upper right-hand portion of the display. The 4D specific parameters are Quality (Q), Volume Angle (A) and Volume Rate (VR). The Status Bar contains instructions on the tasks you can perform at each stage of the 4D imaging process. Remember to take a look at the Status Bar as needed.

4D Touch Panel Overview

Most 4D Touch Panel screens contain some similar controls. Refer to the table below for descriptions of these controls. Controls that are unique to or that contain slightly different functionality are described in the respective sections.

Table 6-34 Common 4D Touch Panel Controls

Preset Parameter	Description
Tile	The display can be divided into 1, 2, or 4 windows for Render view (Render = On) and 1 or 4 windows for Sectional view (Render = Off).
Reset Curve	Resets the three-point curve to a straight line.
Direction	Adjusts the view direction of the ROI.
Visualization	Sectional, Render, VCI, or Tomographic Ultrasound Imaging (TUI). Render view displays one rendered image, or reference image(s) and rendered image.
Focus Position	Adjusts the focal position.
Volume Angle	Sets the range of the volume sweep.
Quality	Balances speed with line density. Max combines the highest density with the slowest speed; Low combines the lowest density with the highest speed.

4D Presets

Real-Time 4D/Static 3D Presets

1. When entering 3D/4D mode, select the **Preset** tab.
2. Select one of the preset settings for data acquisition and display. Presets are defined in the preset file and differ by application.

Table 6-35 Common 4D Touch Panel Controls

Preset Parameter	Description
Over Write	Overwrite the application preset file with the changes you just made.
Create New	Create a new user application preset file based upon the current exam category and application.
User	Used to define new user presets for a given application.

Static 3D Presets

1. When entering 3D/4D mode, select **Static 3D**, then the **Preset** tab.

2. Select one of the preset settings for data acquisition and display. Presets are defined in the preset file and differ by application.

Performing a 4D Scan


Visualizations

4D provides two types of views for displaying and working with images: Sectional, Render, and Tomographic Ultrasound Imaging (TUI).

Sectional View

Sectional view contains one display for each sectional plane.

Figure 6-122 Monitor Display in Sectional View



1. Sectional Image A
2. Sectional Image B
3. Sectional Image C


Render View

The ultrasound system continuously displays the 4D rendered image.

NOTE

When the tile selection is single, only the rendered 4D image appears. When the tile selection is quad, the sectional images are located in 3 quadrants with the rendered 4D image in the fourth.

Figure 6-123 Quad Tile Render View



Orientation Help

When viewing a 4D image on the display monitor, it's sometimes difficult to recognize the orientation. To help, the system displays a three-dimensional drawing to illustrate the orientation. This drawing displays ONLY in sectional view.

Figure 6-124 Orientation Help Graphic



Acquiring and Rendering a 4D VOI

Starting with a 2D Image

To create a 4D image, start with an optimized 2D image. The 2D image serves as the mid-line for the resulting 4D image.

1. Connect the appropriate 4D-compatible probe, leaving the probes in their respective holders.

NOTE

If the appropriate 4D probe is not connected, the original 3D Touch Panel appears.

2. Obtain a 2D image. Optimize the image as usual.

Entering 3D/4D Mode

In 3D/4D mode, the type of scan to perform is chosen: 4D or Static 3D.


1. Press **3D/4D** to enter 3D/4D mode. The first time 3D/4D is pressed, the system is in B Pre mode.

NOTE

The location of the number of focal zones might change when entering 3D/4D mode, since the number of zones is pre-determined by the default ROI.

The default acquisition mode varies by application. In the OB application, the default acquisition mode is Real-Time 4D. For all other applications, Static 3D is the default acquisition mode. When entering pre-mode, an ROI graphic may appear on the monitor display that defines the initial ROI (Region of Interest) of the volume.

Figure 6-125 ROI Graphic



2. Press the **Preset** tab. Select one of the preset settings for data acquisition and display. Presets are defined in the preset file and differ by application.

Acquiring a 4D Volume of Interest (VOI)

To create a 4D image, you start with an optimized 2D image. The 2D image serves as the mid-line for the resulting 4D image.

During 4D image acquisition:

- Frame Averaging is disabled.
- You cannot change the transmit frequency.
- To change the position of the focal zone, adjust the depth of the VOI.
- You cannot change the number of focal zones.

Follow below steps to acquiring a 4D Volume of Interest (VOI).

1. Connect the appropriate 4D probe, leaving the probe in the respective holder.

NOTE

GE HealthCare recommends not to connect 4D straight handle probe to the center probe port. It will prevent user to open the drawer.

2. Select the appropriate 4D probe from the probe indicator.
3. Obtain a 2D image. Optimize the image as usual.
4. Press **3D/4D** to enter 3D/4D mode. An ROI graphic appears. 4D is selected.
5. Define the Volume of Interest (VOI) to be scanned. Press **Pos Size** until **Pos** is highlighted, then move the Trackball to position the VOI. Then press **Pos Size** until the **Size** is highlighted, then move Trackball to re-size the VOI. Only the area defined within the VOI is rendered.
6. Adjust the volume angle. This defines the range of the volume sweep. A small sweep angle results in a lower number of slices with a high volume rate.
7. Set **Render** on to begin 4D acquisition, press **Start** key. You DO NOT have to hold the probe steady during data acquisition.

During data acquisition, you can manipulate the VOI to see different views of the image. To rotate the VOI left or right, use the **PW/Y** control. To rotate the VOI forward or backward, use the **CF/Z** control. To rotate the VOI in a circular motion, use the **M/X** rotary control.

To return to 3D/4D pre-mode, press **Pre Mode**.

NOTE

If the volume size is too large, the message "Volume Size Too Big - Quality Degraded" displays in the status bar. The system changes the quality automatically to below the upper limit and displays the proper Quality value in the information window.

8. To complete the acquisition, press **Freeze** or **Stop**.
9. Store the image.

V-Live 2.0


V-Live 2.0 can be activated when performing static 3D and real-time 4D functions. The lighting source used by the rendering can be positioned by the user around the rendered 3D object on the sphere. By highlighting the side structure, a more realistic three-dimensional image will be obtained.

Image rendering based on subsurface scattering technology is recognized as a mature technology. Most of current surface reconstructions use frontal lighting on rendered objects, which can cause images to look flat. Using the high-resolution volumetric imaging (V-Live 2.0) feature, lighting sources used for rendering can be positioned by the user around the rendered 3D object on a sphere. By highlighting the side structures, the three-dimensional image is significantly improved.

You can activate **V-Live 2.0** during live 4D scanning.

1. While scanning in 4D, select **V-Live 2.0** on touch panel.
2. The flashlight icon displays on the screen. Press **Light** to continue. Move the trackball to adjust the lighting direction for the scanning area.

Figure 6-126 V-Live 2.0 lighting sample




User can also activate **V-Live 2.0** from Render tab after 3D/4D image acquisition.

Advanced Features

1. In **4D Render** or **Static 3D Render** tab, select **V-Live 2.0** from the drop-down menu of **Rend 1 Gray** on touch panel.

Figure 6-127 V-Live 2.0 on Touch Panel



2. The flashlight icon displays on the screen. Press **Light** to continue. Move the trackball to adjust the lighting direction for the scanning area.

NOTE

The lighting area can be adjusted 360 degree.

Sectional VOI Acquisition

Sectional view provides three separate views of the same image: Longitudinal (original image), Transverse (elevational), and Coronal (horizontal).

1. In the 4D tab, Render defaults to On (Render mode). Change Render to Off for Sectional view.
2. To select a reference image, use the Ref Image control on the Touch Panel. The reference image selected contains the focus for control panel keys, allowing you to manipulate or optimize that image.

Table 6-36 4D Data Acquisition Parameters

4D Parameter	Description
Restore View	Resets all parameters back to the original values or chosen presets.
Tile	Selections: Single, Quad. You can divide the display into 1 or 4 windows.

4D Parameter	Description
Visualization	Sectional, Render, or Tomographic Ultrasound Imaging (TUI). Render view displays one rendered image, or reference image(s) and rendered image.
Orientation Help	Displays a three-dimensional drawing to illustrate the orientation. Only displays in sectional view.
Volume Angle	Sets the range of the volume sweep.
Quality	Selections: Max, Hi2, Hi1, Mid2, Mid1, Low. Used to balance speed with line density. Max combines the highest density with the slowest speed. Low combines the lowest density with the highest speed. BQ displays on the display.

Render VOI Acquisition

Rendering allows distinguishing of subtle anatomical detail. All areas of a VOI can be rendered, or just certain regions of the VOI. The region defined for rendering is referred to as the Render Box.

1. Define the area to be rendered. For example, with an image of an entire fetus, only the fetal face may need to be rendered. Therefore, define the fetal face as the VOI.
2. Set Render to On.

Table 6-37 4D Data Acquisition Parameters - Render Mode

4D Parameter	Description
Restore View	Select to reset all parameters back to the original values or chosen presets.
Tile	Selections: Single, Dual, Quad. You can divide the display into 1, 2, or 4 windows.
Visualization	Sectional, Render, or Tomographic Ultrasound Imaging (TUI). Render view displays one rendered image, or reference image(s) and rendered image.
3D Orient	When selected, changes the orientation of the image on the monitor display. Selections include: 0 degrees, 90 degrees, 180 degrees, and 270 degrees.
Volume Angle	Sets the range of the volume sweep.
Quality	Selections: Max, Hi2, Hi1, Mid2, Mid1, Low. Used to balance speed with line density. Max combines the highest density with the slowest speed. Low combines the lowest density with the highest speed.
Activate Curve	Define a three-point curved surface for the render window using the Trackball.
Reset Curve	Reset the three-point curve to a straight line.
Mix	Selections: 0-100% in increments of 2. Allows you to mix a Rend Mode 1 mode with a Rend Mode 2 mode. Always select two modes.

4D Parameter	Description
Threshold	Selections: 0-255. Sets a lower threshold below which weaker echoes are removed.

3. Select the **Render Setting** tab.

The Render Setting tab allows selection and combining of gray-scale and color rendering modes.



HINTS

When using Surface modes, adjust the Lower Threshold to recognize border structures more clearly.

Table 6-38 4D (Data Acquisition) Render Parameters

4D Parameter	Description
Direction	The ROI determines the region that is rendered during 4D acquisition. The direction in which this ROI is viewed can be changed. Selections: Up/Down, Down Up, Left/Right, Right/Left, Front/Back, Back/Front.
Gray Map	Displays the gray map selections on the display monitor. Select maps using the Trackball.
Colorize	Displays the tint map selections on the display monitor. Select maps using the Trackball.
Render Mode	Select Gray or Inversion. Selecting Inversion inverts the gray values of the rendered image (e.g., image information that was black becomes white, and vice versa).

4D Parameter	Description
Render 1 / Render 2	Allows you to combine render mode values from render mode 1 and render mode 2. Select the render map combination from the upper-left portion of the monitor display. Select map combinations using the Trackball. Render Mode 1 Selections: Surface Smooth, Surface Texture, Transp Max, Transp X-ray, Transp Min. Surface Smooth - Surface displays in a smoothed texture mode, which means that the gray values of the surface are identical with the gray values of the original 2D scan. Surface Texture - Surface displays in texture mode, which means that the gray values of the surface are identical with the gray values of the original 2D scan. Transp Max. - Displays the maximum intensity of gray values in the ROI. This is helpful for viewing bony structures. Transp X-Ray - Displays all gray values in the ROI. Trans Min. - Displays the minimum number of gray values in the ROI. This is helpful for viewing vessels and hollow structures. Render Mode 2 Selections: Surface Smooth, Light, Gradient Light, Transp Max, Transp X-ray, Transp. Min. Surface Smooth - Surface displays in a smoothed texture mode, which means that the gray values of the surface are identical with the gray values of the original 2D scan. Light - Surface displays in light mode. Structures in the near field are brighter; structures in the far field are darker. Gradient Light - Surface displays as if it is illuminated from a spot light source. This is helpful if the displayed surface is surrounded by hypoechoic structures (for example, liquids). Transp Max. - Displays the maximum intensity of gray values in the ROI. This is helpful for viewing bony structures. Transp X-ray - Displays all gray values in the ROI. Transp Min. - Displays the minimum number of gray values in the ROI. This is helpful for viewing vessels and hollow structures.
Adv. Rend	Surface Smooth Plus - Surface displays in more smoothed mode. Transparence Plus - Makes the image/information more transparent.
Transparenc y	Selections: 20 to 250. Sets the transparency of the image. The higher the number, the more transparent the gray scale information.

Manipulating the Volume of Interest (VOI)

The 3D/4D ROI is a tangible anatomical object that can be viewed and manipulated easily using the **Trackball** and other control panel keys.

If the monitor display is in Sectional view, select the desired reference image before manipulating the image.

Moving through the VOI

To move through the image to view a particular slice, press **Shift**.

This allows a displacement of the center of rotation along the intersection lines of the sectional planes A, B, and C. The displacement of the center of rotation leads to the display of parallel sectional images.

Zooming the Image

Rotate **Zoom** to zoom on the image.

Moving the VOI Position

To move the position of the VOI, press **Scan Area** button on the control panel until **Pos** is highlighted, then move the **Trackball** left, right, up and down, as needed.

Resizing the VOI

To re-size the VOI, press **Top Right Key** on the control panel until **Size** is highlighted, then move the **Trackball**.

Stopping 4D Image Acquisition

To stop acquiring a 4D image, press **Freeze/Stop**.

4D VOI Post-Processing

When you press **Freeze**, the system menu may differ, depending on whether you are in Render view or Sectional view.

Volume CINE

The system constantly stores CINE images to play back and review those images. CINE is useful for focusing on images during the specific part of the heart cycle or to view short segments of a scan session.

To activate CINE in 4D:

1. Press **Freeze**.
2. Select the **VolCine** tab.

Table 6-39 4D Cine Parameters

Preset Parameter	Description
Loop Mode	Selections include: One Way, BiDirectional (two-way). One Way - plays one loop sequence forward. BiDirectional - plays the sequence forward and backward.
First	Displays the first volume in the CINE loop.
Last	Displays the last volume in the CINE loop.
Run/Stop	Starts and stops the CINE loop.
Loop Speed	Adjusts the CINE loop speed.
Volume by Volume	Used to select an individual volume in the CINE loop.

3. If you were in Render visualization mode when you entered 4D CINE mode, you can also press **Pre Mode** Key to return the Pre Mode.

If you were in Sectional visualization mode when you entered 4D CINE mode, you can press **Pre Mode** key to return the Pre Mode.

4. To re-start real-time 4D acquisition, press **Freeze**.

Static 3D

A single sweep, single volume static 3D image can be created.

Performing a Static 3D Scan

1. Connect the appropriate 4D probe, leaving the probe in the respective holder.
NOTE
GE HealthCare recommend not to connect 4D straight handle probe to the center probe port. It will prevent user to open the drawer.
2. Select the appropriate 4D probe from the probe indicator.
3. Obtain a 2D image. Optimize the image as usual.
4. Press **3D/4D**.
5. Press Static 3D. Set Render to On.
6. Define the Volume of Interest (VOI) to be scanned. Press **Pos Size** until **Pos** is highlighted, then move the Trackball to position the VOI. Then press **Pos Size** until the **Size** is highlighted, then move **Trackball** to re-size the VOI. Only the area defined within the VOI is rendered.
7. Adjust the volume angle. This defines the range of the volume sweep. A small sweep angle results in a lower number of slices with a high volume rate.
8. Make sure the probe is held steady. Press the **Start** key to start acquisition.

NOTE

During 3D acquisition, the CW, PW, M, and Depth control panel keys are not available.

NOTE

When the 3D acquisition begins, the Touch Panel appears blank for a brief moment.

9. Hold the probe steady until the system stops automatically. You will know the acquisition has stopped when the system menu changes to display the Render Setting, 3D Rotational Cine, and Scalpel tabs.
10. Store the image.

You can also activate **V-Live 2.0** to use the light during Static 3D Scan. Refer to *V-Live 2.0* on page 520 for detail information.

Static 3D Sectional View

Table 6-40 3D After Acquisition Parameters - Sectional View

Preset Parameter	Description
Orientation Help	Displays a 3-dimensional drawing to illustrate the orientation. Only displays in sectional view.

Static 3D Render View

Table 6-41 3D After Acquisition Parameters - Render View - Page 1

Preset Parameter	Description
Edit/Accept ROI	Selections include Edit, Accept. Edit - Select to adjust the size of the Region of Interest (ROI). Accept - accepts the active 3D image.

Preset Parameter	Description
3D Orient	When selected, changes the orientation of the image on the monitor display. Selections include: 0 degrees, 90 degrees, 180 degrees, and 270 degrees.

Scalpel

Scalpel allows editing/cutting sections of a 3D image. Scalpel is available only on a rendered image.

1. Select **Scalpel**.

Table 6-42 Scalpel Parameters

Preset Parameter	Description
Cut Mode	Selections: Inside Contour, Outside Contour, Inside Box, Outside Box, Eraser Big, Eraser Small. Inside Contour, Outside Contour - Allows tracing the portion of the image to cut. Trace Outside removes all portions of the image that fall outside the traced region. Trace Inside removes all portions of the image that fall inside the traced region. Inside Box, Outside Box - Displays a box used to define the portion of the image to cut. Outside Box removes all portions of the image that fall outside the box. Inside Box removes all portions of the image that fall inside the box. Eraser Big, Eraser Small - Provides a big and small eraser tool used to define the portion of the image to cut by hand. Available only if Depth is Full.
Cut Depth	Selections: Full, Define. Full -The entire depth of the selected region will be cut. Define - Allows you to define the depth to cut using the Depth control panel knob.
Undo Last	Undoes the last cut only.
Redo	Select to redo scalpel.
Undo All	Undoes all cuts since you entered Scalpel mode.
Done	Applies to User Defined Cut Depth when complete.

2. Select the cut mode.
3. Use the **Trackball** and Set key to define the portion of the image to cut. Press Set to start, move the **Trackball** to define the region, then press Set again to cut the image. The portion is removed.

To undo the last cut, select **Undo Last**.

To undo all cuts in the current session, select **Undo All**.

NOTE

With the cut image displayed, if you attempt to switch to the Static 3D tab to edit the ROI, the following warning message appears: Scalpel changes will be lost. Do you want to continue? [Yes/No].

3D Rotation CINE

3D Rotation CINE allows you to view the 3D image from various angles.

To activate rotation CINE in 3D:

1. Press **Freeze**.
2. Select the 3DRot Cine tab.

Table 6-43 3D Rotation Cine Parameters

Preset Parameter	Description
Rotational Angle	Sets the rotational angle of the 3D image over which the CINE loop is played. Typical values are 30, 45, 60, 90, 180 and 360 degrees.
Step Angle	Sets the step angle between individual frames in the CINE loop.
Rotation Axis	Sets the axis about which the CINE loop is calculated. Selections X and Y.
Loop Mode	Selections include: One Way, BiDirectional (two-way). One Way - plays one loop sequence forward. BiDirectional - plays the sequence forward and backward.
First	Displays the first volume in the CINE.
Last	Displays the last volume in the CINE.
Run/Stop	Starts and stops the CINE sequence.
Start Angle	Used to select the starting angle in the CINE loop range. The default Start Image is calculated from the rotational angle as: $-1 \times \text{Rotational angle} / 2$. If you adjust the Start Image, the Rotational Angle is re-set to be the value of the adjusted Start Image.
End Angle	Used to select the ending angle in the CINE loop range. The default End Image is calculated from the rotational angle as: $\text{Rotational angle} / 2$. If you adjust the End Image, the Rotational Angle is re-set to be the value of the adjusted End Image.
Image by Image	Used to select an individual image in the CINE loop.

VOCAL

Use VOCAL (Virtual Organ Computer-aided Analysis) to visualize and calculate the volume of anatomical structures, such as a tumor lesion, cysts, and the prostate. VOCAL is available after a Static 3D or Real-Time 4D acquisition.

1. Select **Vocal**. Specify the volume calculation method (**Manual**, **Contour Detect**, **SemiAuto Detect**, or **Sphere**). Select the reference image you want to use to perform the trace by selecting **Ref Image A**, **B**, or **C**. Press **Start**.
2. Trace the anatomy using the **Trackball**. Press **Set** to start and end the trace. You must go across the dotted line for the trace to take effect (it turns yellow). The trace is performed on each image slice, separated by the rotational step angle. Rotate the **Rot. Ref** dial until you have completed the total of the required rotations (for example, if you've selected 30 degrees, you need to complete six traces if you've selected Manual). After you've completed the trace target, the **Calc Volume** button is active for you to press. The calculated VOCAL image appears in the lower, right-hand corner of the display. You are now in the edit state.

NOTE

Trace not used for Sphere. For Sphere, set the Poles.

3. Edit as necessary. You can apply a shell, adjust its thickness, navigate through the reference angles, or restart the VOCAL.

VOCAL Touch Panel states are described in the following tables.

Table 6-44 VOCAL Parameters on Setup Touch Panel

Parameter	Description
Manual	When you select the Manual method, you need to perform a manual trace on each of the rotation angles.
Contour Detect	When you select the Contour Detect method, you need to perform a manual trace on each of the rotation angles.
SemiAuto Detect	When you select the Semi Automatic Detect method, you need to perform a trace on only two rotation angles. The system applies an algorithm to define the traces.
Structure	Structure is only available with the SemiAuto Detect method. Select Hypo, Cystic, or Hyper/Iso.
Sphere	Calculates the volume based on the pole settings.
Pole 1	Adjust the upper contour point (green arrow) of the structure.
Pole 2	Adjust the lower contour point (green arrow) of the structure.
Rotational Step Angle	Specify the angular spacing between contour traces. Typical values are 6, 9, 15, and 30 degrees. The number of planes varies by this formula: $180 \text{ degrees} / \text{selected rotational step angle}$.
Ref Image	Use this to select the image you want to use to perform the trace.
Start	Press Start when you're ready to perform the trace.
Rot.Ref ### Back/Next	Select Next/Back to move to the next image for contour definition in the rotation step.

Table 6-45 VOCAL Calculate Volume Touch Panel

Parameter	Description
Calc Volume	Press Calc Volume to initiate the VOCAL image calculation.
Clear	Press Clear to remove the trace from the image.
Restart Vocal	Press Restart Vocal to return to the initial VOCAL state.

The Shell Modes allows you to construct a shell or contour “around” the structure of interest, which enables you to distinguish between the contour of the targeted structure and the contours the inside and outside of the structure.

Table 6-46 VOCAL Parameters on Edit Touch Panel


Parameter	Description
Shell Off	Select Shell Off if you do not want a shell around the VOCAL image.
Inside	Select Inside if you want a shell inside the volume.
Outside	Select Outside if you want a shell outside the volume.
Symmetric	Select Symmetric if you want half of the shell thickness inside and half outside the volume's perimeter.
Shell Thickness	Adjust to vary the thickness of the shell.

Tomographic Ultrasound Imaging (TUI)

Tomographic Ultrasound Imaging (TUI) is a visualization mode which presents data as parallel slices (planes) through the dataset. This method of visualization is consistent with CT and MRI. The distance between the different planes can be adjusted.

1. Select **TUI** as the Visualization mode.
2. Start acquisition.
3. If in 4D, press **R** to end the acquisition. This step is not required in Static 3D.
The reference image and the number of specified slices appears. The reference image always displays and indicates which slices you are currently viewing as solid lines.

Figure 6-128 TUI 3x3 Example



a. The TUI reference image that shows the slice position. This image is orthogonal to the reference image.

- b. The number and green asterisk shows the position of each slice. A green asterisk indicates the center image (A, B or C-plane).
- c. Slice distance displays when the slices are in certain intervals.
- d. A solid line indicates the slice appears on the monitor.

A dotted line indicates the slice did not appear on the monitor.

NOTE

TUI with Color is only available with Static 3D, not with 4D.

- 4. Adjust the number of slices and slice distance.

You can adjust the number of slices by using the **Slices** rotary. You can adjust the distance between the slices using the **Slice Distance** rotary. Max value is 40mm.

- Select Display Format from 1x1, 1x2, 2x2 and 3x3.
- Change the center image via **Ref. Image** if needed (Reference image A, B or C).
- Move forward/backward through the slices via Prev./Next Slice.
- The following features are supported in TUI: Rotation (X/Y/Z), Trackball (Move the position) and Gain.

Whizz RenderLive

Whizz RenderLive helps to find the render start position to easily separate solid tissue in front of the rendered object.

The Whizz RenderLive algorithm “looks” for the transition from solid tissue to liquids and positions the “Render Start” into the liquid area visualized by the green render start line. The render start line is not a straight line but a “free” trace for optimal adaptation to the rendered object.

NOTE

The render start line can be positioned by pressing auto. It is possible to adjust the line trace manually but not its sensitivity. For more flexibility/sensitivity, Whizz RenderLive has to be activated.

Using Whizz RenderLive

1. Start the Render Visualization Mode.
2. Press the **Whizz RenderLive** knob.
3. To adjust the distance between the render start position and the render object, rotate the **Sensitivity** control below the touch panel. A high value indicates a smaller distance.

NOTE


In the case that Whizz RenderLive is not used, the Render Start line can also be modified manually. Press the trackball button Curv to activate Curved Render Start and move the trackball to modify the line.

Advanced Features

NOTE

Whizz RenderLive is not available in Sectional or TUI Mode.

Figure 6-129 Whizz RenderLive example



Whizz Follicle

Whizz Follicle helps detect low echogenic objects (e.g. follicles) in an organ (e.g. ovary) and analyzes their shape and volume. From the calculated volume of the object, an average diameter is calculated. All objects detected using whizz follicle will be listed according to size. The calculation results are displayed in the bottom left corner, with the object numbers color-coded to denote the object on the image. If the mouse cursor hovers over a specific item on the list, the respective object in the image is highlighted (and vice versa). The color of the object is bound to its position on the list.

On the acquired three-dimensional image data, low-echo object areas (such as follicles) in organs (such as ovaries) are detected by manual cursor selection or manual measurement of height and length, and the trace and number will be marked on the screen. The system

Advanced Features

calculates the volume and average diameter of the object based on the measurement results, and display them in the Whizz Follicle table from largest to smallest based on the results.

Figure 6-130 Whizz Follicle Display Screen




Figure 6-131 Whizz Follicle Touch Panel




Table 6-47 Whizz Follicle Parameters on Touch Panel

Parameter	Description
SemiAuto	Each object is selected manually with the system cursor but traced / measured automatically.
Manual	Each object has to be selected and measured manually.

Parameter	Description
Angular	The box has the shape of a rectangle.
Rounded	The box has an elliptic shape with rounded corners.
Left Ovary	Choose to measure Left Ovary with GYN application.
Right Ovary	Choose to measure Right Ovary with GYN application.

Using Whizz Follicle

Manual

1. Create a 3D Static Volume of the desired organ.
2. Press the **Volume Analysis** in **Static3D Render** tab.
3. Select **Whizz Follicle** on the Touch Panel
4. Select **Manual** and adjust the ROI shape, if desired.
5. Start the measurement by selecting either **Left Ovary** or **Right Ovary**.
6. The green measurement cross appears on the monitor.
7. Perform measurements:
 - a. Position the start-point of the long diameter with the Trackball and press **Set**.
 - b. Position end-point of the long diameter with the Trackball and press **Set**.
 - c. Position the start-point of the short diameter with the Trackball and press **Set**.
 - d. Position end-point of the short diameter with the Trackball and press **Set**. The second measurement displays.
 - e. To start the next measurement move the Trackball and repeat step 7.
8. Press **AddToReport** on the Touch Panel.

NOTE
A message will appear on the screen to indicate the report is added successfully.
9. User can press **New Analysis** and start a new analysis by repeating step 5 to step 8.

Semi Auto

1. Create a 3D Static Volume of the desired organ.
2. Press the **Volume Analysis** in **Static3D Render** tab.
3. Select **Whizz Follicle** on the Touch Panel
4. Select **SemiAuto** and adjust the ROI shape, if desired.
5. Start the measurement by selecting either **Left Ovary** or **Right Ovary**.

6. Move the green arrow cursor to the middle of the desired area to measure and press **Add/Rem.** trackball key. Then use trackball keys to add, remove, cut or merge the desired area.

Figure 6-132 Whizz Follicle Trackball keys



- | | |
|-----------------|--|
| Add/Rem. | Position the system cursor over the detected object, then press Add/Rem. key to add or remove objects manually. |
| Cut | Position the system cursor over the detected object, then press Cut key and use the trackball to trace the area you want to remove and press Cut key again to cut the traced area. |
| Merge | Position the system cursor over the detected object and press Merge key and move the trackball to another object. Then press Merge key again to combine these 2 objects. |

7. Press **AddToReport** on the Touch Panel.
8. User can press **New Analysis** and start a new analysis by repeating step 4 to step 7.

User can also use the **Add Manual** key while **SemiAuto** is selected when objects are not detected or wrongly measured.


1. Press **Add Manual** button.
2. The measurement cross appears.
3. Perform measurements:
 - a. Position the start-point of the long diameter with the Trackball and press **Set**.
 - b. Position end-point of the long diameter with the Trackball and press **Set**.
 - c. Position the start-point of the short diameter with the Trackball and press **Set**.
 - d. Position end-point of the short diameter with the Trackball and press **Set**. The second measurement displays.
 - e. To start the next measurement move the Trackball and repeat step 3.

Whizz Follicle Report

1. Press **Report** on touch panel.
2. Select **Whizz Follicle** in **Worksheet** Tab.

3. Press **Whizz Report** to generate report.

Figure 6-133 Whizz Follicle Worksheet and Report Touch Panel



Whizz Note

Whizz Note allows the system to display folders and files from removable media. The user can import and display useful clinical datasheet for reference during scanning.


1. Insert a removable media to the ultrasound system.

NOTE

Supported media includes USB, HDD, CD/DVD and Network Storage.

2. Press **Assistant** key on control panel.
3. Click **Whizz Note** on the touch panel.

Figure 6-134 Note on Touch Panel




4. A window displays on the screen. Select desired file, click **Import** to import files to **Whizz Note**.

NOTE

Whizz Note only supports importing files in bmp, png, jpg, jpeg, gif, ics or pdf format. File sizes should be less than 2 Megabytes each.


Figure 6-135 Import files to Whizz Note



Advanced Features

- The imported files will be listed in the **Whizz Note** window. User can select desired files and click **Show** to open the file as a reference during operating the system. You can also select **+** icon to import more files.

Figure 6-136 Show the file in Whizz Note



- User can manipulate Whizz Note window while operating the system.

Figure 6-137 Whizz Note Main page




Table 6-48 Controls on Whizz Note main page

No	Control	Description
1	Page control	Select to display previous or next page.
2	Minimize/Maximize	Select to maximize or minimize the window.
3	Close	Select to close Whizz Note.
4	Size	Move the cursor till it turns to an arrow, hold and drag to change the size of the window.

7. To add more files, click the + icon on the right upper corner.

VCO for customer

VCO for Customer allows a remote expert(s) who has (have) been authorized by the user to connect to ultrasound system remotely. It enables remote imaging support, communication and application training through the facility server and internet connection between the 2 components of the feature, DClient (Device Client) and EClient (Expert Client). DClient is integrated in the software of the ultrasound device. EClient is a software application that can be installed on a personal computer or mobile device.



WARNING

VCO for Customer is not intended for diagnostic use since GE HealthCare cannot validate the shared image feed on the EClient side.



CAUTION

VCO for Customer is not available during biopsy procedure or cardiac application.

Outside Hospital version

Before using VCO for Customer (Outside Hospital version)

Before using VCO for Customer, please follow below steps to prepare both DClient and EClient sides.

1. Make sure internet configuration for ultrasound system is completed before using VCO for Customer.
2. Please connect USB2.0 drive-free camera and microphone through USB ports on ultrasound system to share with video and audio during VCO for Customer session.
3. Enter <https://jysx.gehealthcloud.cn> into internet browser of local device with assigned admin account information to create user account(s) for EClient expert(s).

NOTE

Admin account assigned by GE HealthCare can be used to create additional ultrasound system (DClient) accounts as well as EClient accounts. Each expert needs a separate account.

4. Share created account information (username and password) with each EClient expert(s).
5. Inform EClient expert(s) to download and install EClient application on the device they plan to use during the sharing session.
 - For Windows devices, use internet browser to enter <https://jysx.gehealthcloud.cn> for downloading and installing EClient application.
 - For Android mobile devices, scan below barcode to download and install EClient application.




- For iOS mobile devices, visit App Store and search for 技影随行-远影 to download and install EClient application.
6. You can start using VCO for Customer now.

Using VCO for Customer (Outside Hospital version) - DClient Side

1. Press **Assistant** key on the control panel, select **VCO for Customer** on the Touch Panel.

Figure 6-138 Enter VCO For Customer



2. Select **Outside Hospital Version** and **Region** in the pop up window, then press **Confirm** to continue.

Figure 6-139 Select Outside Hospital version



NOTE

Please select appropriate region from the pop up window for the first time login the **Outside Hospital Version**.

3. Input username and password on login page.

NOTE

Only authorized user can use VCO for Customer.

NOTE

User name and password can be assigned by the hospital administrator. If password is lost, please contact administrator or GE HealthCare service to reset password.

NOTE

If connection error occurs, please check internet connection first. If internet connection is good, please check if server address is correct.

4. Press **Login** to enter VCO for Customer. User can select desired scenario to start real-time assistant, training education, quality control and/or offline communication.

Figure 6-140 Outside Hospital version scenario



NOTE

When server connection is not good, there will be a caution message says server connection is poor.

Below table lists the description for the scenario on home page.

Table 6-49 VCO for Customer home page on ultrasound system

Scenario	Description
Real-time Assistant	• DClient user requests to the remote expert to help provide support/guidance during imaging. • DClient user shares the ultrasound scan screen and can have a live conversation. The video camera can be placed in position to allow the EClient user to see imaging actions like probe position of the DClient user. • EClient user guides user of DClient with imaging operation through live discussion.

Scenario	Description
Quality Control	• DClient user enters quality control mode. • EClient user(quality control center) selects device(s) for the quality assessment to be performed on. • EClient user views multiple scan screens of clients device simultaneously. • EClient user can request live discussion during the quality assessment.
Training Education	• DClient user can invite EClient to training session. • DClient user shares ultrasound scan screen and positions external video camera to allow EClient to view live scanning. • DClient user and EClient users can communicate live through Video and Audio.
Offline Communication	• DClient user sends JPG or wmv images to EClient for communication. • EClient user can view images and send text message between expert clients.
GE Online	• DClient user and EClient user can test the connection or start online training. • DClient user can consult application issue during using the feature.

Real-time Assistant

To use Real-time Assistant:

1. Select **Real-time Assistant** on home page.



WARNING

You are about to activate remote sharing function. Please pay attention to the protection of patient information to avoid leakage

2. Select expert(s) from **Experts List** and press **Enter Room** to start Real-time Assistant.

NOTE
Experts list is generated automatically after user account(s) are created in step 4 of section "Before using VCO for Customer".

Figure 6-141 Experts List



CAUTION
Please stop sharing when VCO is running at non-scanning window to prevent any potential PHI (Patient Health Information) disclosure.

Main user interface controls in Real-time Assistant scenario are illustrated as below:



Table 6-50 User Interface controls in Real-time Assistant

Icon No.	Icon Name	Description
1	Audio	Enable/disable audio output.
2	Video	Enable/disable video output.
3	Share	Enable/disable screen sharing.
4	Setting	Change the setting for image and other configuration.
5	Exit	End meeting and exit VCO for customer.
6	Hide	Minimize the main control menu to the top of the screen.
7	Local	Video Window on DClient side.
8	Remote	Video Window on EClient side.
9	Members	Online members list. User can mute all other users on the list in this window.
10	Sharing area	Screen sharing area.

Training Education

DClient user can invite EClient user to training session or share ultrasound scan screen and position external video camera to allow EClient user to view live scanning.

DClient user and EClient expert can communicate live through Video and Audio. The remote training classroom can support up to 300 users online at the same time.

1. Select **Training Education** on home page.

2. Select expert(s) from **Experts List** and press **Enter Room** to start Training Education.


Main user interface in Training Education scenario is the same as Real-time Assistant, refer to *Real-time Assistant* on page 548 for detail information.

Quality Control

Select **Quality Control** on home page to enter Quality Control mode.

An expert can see and communicate with multiple devices in Quality Control mode. For the detail information of user interface in Quality Control scenario, please refer to *Real-time Assistant* on page 548 for more information.

Figure 6-142 Quality Control key



User can also press the key a on the quality control page to request Quality Control from EClient expert.

Offline Communication

When the EClient expert is not available for real-time assistant, user can send images/cine loops for offline communication.

Follow below steps for offline communication:

1. Complete the exam and exit VCO for Customer.
2. For traditional workflow, select the exam and click **Active Image**. For simplified workflow, select the exam and click **Image Review**.
3. For traditional workflow, select the images/cine loops needs to be sent, then click **'SaveAs' Images**. For simplified workflow, select the images/cine loops needs to be sent, then click **Send To**.
4. For traditional workflow, select **Transfer to VCO For Customer** under the drop-down menu of **Save in archive**. For simplified workflow, select **VCO For Customer** in the pop-up window.

Advanced Features

- For traditional workflow, click **Save**. For simplified workflow, click **Send**. Then the system will transfer the selected images/cine loops without PHI and turn on VCO for Customer automatically.



CAUTION

The patient information should be hidden when using VCO for Customer to transfer images from DClient (Device Client) to EClient (Expert Client).

Figure 6-143 Send Image/Cine loops under Traditional Workflow





Figure 6-144 Send Image/Cine loops under Simplified Workflow



- Input username and password to login VCO for Customer.

7. Select expert from **Experts List(a)** and select files need to sent under **Files(b)**, then click **Send** to send files offline.

NOTE

User can send maximum 5 files in each offline communication message.

Figure 6-145 Send Offline Communication



8. Click **X** on the right upper corner to exit offline communication.

Inside Hospital version

Before using VCO for Customer (Inside Hospital version)

Before using VCO for Customer, please follow below steps to prepare both DClient and EClient sides.

NOTE

Please make sure connecting to the intranet for the system internet configuration.

1. Make sure internet configuration for ultrasound system is completed before using VCO for Customer.


Advanced Features

2. Please connect USB2.0 drive-free camera and microphone through USB ports on ultrasound system to share with video and audio during VCO for Customer session.
3. Configure **VCO for Customer** to an ultrasound system User-Defined key in **Utility > User Configurable Key > User Defined Key**.
4. You can start using VCO for Customer now.

Using VCO for Customer (Inside Hospital Version) - DClient Side

1. Press **Assistant** key on the control panel, select **VCO for Customer** on the Touch Panel.

Figure 6-146 Enter VCO For Customer



2. Select **Inside Hospital Version**, press **Confirm** to continue.

Figure 6-147 Select Inside Hospital version



3. Input username and password on login page.

NOTE

Only authorized user can use VCO for Customer.

NOTE

User name and password can be assigned by the hospital administrator. If password is lost, please contact administrator to reset password.

NOTE

If connection error occurs, please check internet connection first. If internet connection is good, please check if server address is correct.

4. Press **Login** to enter VCO for Customer.


Below table lists the description for the scenario on home page.

Table 6-51 VCO for Customer (Inside Hospital Version) home page on ultrasound system

Scenario	Description
Real-time Assistant	• DClient user requests to the remote expert to help provide support/guidance during imaging. • DClient user shares the ultrasound scan screen and can have a live conversation. The video camera can be placed in position to allow the EClient user to see imaging actions like probe position of the DClient user. • EClient user guides user of DClient with imaging operation through live discussion.

5. Select **Real-time Assistant** to continue.

Figure 6-148 Inside Hospital version scenario

**WARNING**


You are about to activate remote sharing function. Please pay attention to the protection of patient information to avoid leakage

NOTE

When server connection is not good, there will be a caution message says server connection is poor.

- Select expert(s) from **Experts List**, the selected experts will be displayed on the right column. Press **Confirm** to enter the device connecting page.

Figure 6-149 Select Experts



CAUTION

Please stop sharing when VCO is running at non-scanning window to prevent any potential PHI (Patient Health Information) disclosure.

- On the device connection page, make sure that the camera, microphone, speaker, and network equipment are connected, click the **Enter Session** button to enter the meeting. User can set the camera and microphone to be turned on/off.
- Main user interface controls in Real-time Assistant scenario of Inside Hospital Version are illustrated as below:

Figure 6-150 User Interface for Inside Hospital Version




Table 6-52 User Interface controls in Real-time Assistant

Icon No.	Icon Name	Description
1	Microphone	Enable/disable audio output.
2	Video	Enable/disable video output.
3	Screen Share	Enable/disable screen sharing.
4	Setting	Change the setting for image and other configuration.
5	End Conversation	Exit real-time assistant.
6	Close Application	Exit the application.
7	Hide	Hide the menu.
8	Local	Video Window on DClient side.
9	Remote	Video Window on EClient side.
10	Members	Online members list.
11	Sharing area	Screen sharing area.

Using VCO for Customer - EClient Side

Before using VCO for Customer on EClient side, scan the barcode below to get the supplementary instruction.

Figure 6-151 Barcode for downloading supplementary instruction



For detail information on using VCO for Customer on EClient side, please read and follow the supplementary instruction downloaded.

Digital Expert Connect

Digital Expert Connect is intended to facilitate remote non-diagnostic image viewing and review, consultation, guidance, support, and education in real time between ultrasound users and remote users. This device is intended for use by healthcare professionals in a clinical support setting on compatible workstations and mobile devices.

NOTE

This tool is not for clinical diagnostics purposes.

Refer to the Digital Expert Connect User Manual for information on setting up and using Digital Expert Connect.

NOTE

Digital Expert Connect is not available in all regions.

Intensity Ratio

Intensity Ratio(IR) is a generic measurement to calculate the ratio of two intensity (Echo Level, displays EL_dB) measurements. The formula is, Intensity Ratio = Intensity A/ Intensity B.

Purpose

Intensity Ratio is intended to calculate the intensity ratio of two sample areas on a frozen B-mode image.

Probes and presets

Intensity Ratio is available through specific probes and presets. Refer to *Table 5-8 Probe modes of operation and Features* on page 330 for detail information.

Scan Mode

Intensity Ratio is only available on 2D B-mode.

Activation

1. Obtain a frozen B-mode image.
2. Press **Measure** on control panel, select a category on Touch Panel.
3. Select **Cardiac Generic** in Cardiology category or **Generic** in other categories, then press **Intensity Ratio** to enter the page.


IR Calculation

IR can be activated from generic measurement menu on a frozen image. Sample area can be drawn freehand by trace or circle. User shall select the format before using Intensity Ratio.

Advanced Features


Select **Double Trace** or **Double Circle** on **Intensity Ratio** main page directly.

Figure 6-152 Select ROI format on Touch Panel



Draw the sample area by trace if **Double Trace** is selected.

Figure 6-153 ROIs drawn by trace



Draw the sample area by circle if **Double Circle** is selected.

Follow below steps to adjust the ROI for Double Circle :


1. Rotate **Set Sample Area Shape** knob.

Figure 6-154 Intensity Ratio Sub Menu




2. Adjust the **Height**, **Width** and **Tilt angle** for **Double Circle** ROI. Click **OK** to start drawing.

Figure 6-155 Adjust Double Circle ROI



Use trackball and **Set** key to draw by **Double Circle**.

Figure 6-156 ROIs drawn by circle



Display Intensity **A** and **B** values and **IR** in measurement window after performing.


IR Calculation with Follow Up Tool

IR measurement can be activated/calculated when entering Follow Up tool.

1. Select **Auto Copy QA Measurement** from **Utility > System > System Imaging > Follow Up Tool**.
2. Obtain two ROIs, save it as Raw DICOM image, press **Follow Up** button, select the RawDICOM image, click **Freeze** button, press **Measure**, then the samples area (two

ROIs) of intensity 1 and intensity 2 shall be copied to side by side Raw DICOM image comparison from a patient's previous exam to the patient's current active exam.

Figure 6-157 IR Calculation with Follow Up Tool



Worksheet and Report

IR calculation can be transferred to worksheet and report. Press **Report > Whizz Report** to display the IR calculations and report.

NOTE

Measurements taken during the active exam are transferred to the current exam's Worksheet only; previous exam measurement information is not included on active worksheets.

US LI-RADS (ACR)

The Ultrasound Liver Imaging Reporting And Data System (US LI-RADS® v2017) is:

- A standardized system for imaging technique, interpretation, reporting, and data collection for screening or surveillance ultrasound exams in patients at risk for developing HCC.
- Supported and endorsed by the American College of Radiology (ACR).
- Developed by a consortium of diagnostic radiologists and hepatologists with expertise in hepatobiliary ultrasound, with input from and approval by the LI-RADS Steering Committee.

US LI-RADS US requires two types of assessments, both pertaining to the whole liver rather than specific nodules or other observations:

- US category
- US visualization score

US Category

The US Category summarizes the main results and helps determine the most appropriate follow-up. These scores are possible:

- US-1 Negative
- US-2 Subthreshold
- US-3 Positive

Table 6-53 US Category Details

Category	Concept	Definition
US-1 Negative	No US evidence of HCC	No observation OR Only definitely benign observation(s)
US-2 Subthreshold	Observation(s) detected that may warrant short-term US surveillance	Observation(s) < 10mm in diameter, not definitely benign
US-3 Positive	Observation(s) detected that may warrant multiphase contrast-enhanced imaging	Observation(s) ≥ 10mm diameter, not definitely benign OR New thrombus in vein

US Visualization Score

The **US Visualization Score** reflects technical or other factors that may affect liver visualization or nodule detection. This information helps to communicate the expected level of sensitivity of the screening exam for HCC detection in an individual patient. Data on visualization scores may be used for quality assurance and to inform future refinements of LI-RADS and LI-RADS-related management guidelines. Three visualization scores are possible:

- A. No or minimal limitations
- B. Moderate limitations
- C. Severe limitations

Table 6-54 US Visualization Score Details

Score	Concept	Examples
A. No or minimal limitations	Limitations if any are unlikely to meaningfully affect sensitivity	Liver homogeneous or minimally heterogeneous. Minimal beam attenuation or shadowing. Liver visualized in near entirety.
B. Moderate limitations	Limitations may obscure small masses	Liver moderately heterogeneous Moderate beam attenuation or shadowing Some portions of liver or diaphragm not visualized
C. Severe limitations	Limitations significantly lower sensitivity for focal liver lesions	Liver severely heterogeneous Severe beam attenuation or shadowing Majority (>50%) of liver not visualized Majority (>50%) of diaphragm not visualized

User Requirement

Probes and Presets

LI-RADS(ACR) is only offered in Abdomen application using ABD or ABD_PEN presets.


Activation

There are two entrance to activate LI-RADS(ACR). Operator can use Measure button to activate:

Advanced Features


1. Press **Measure**, select **Abdomen** category, press **LI-RADS(ACR)** on the Touch Panel to activate LI-RADS(ACR).

Figure 6-158 LI-RADS on Touch Panel



2. Select **US category** and **US visualization score**, input comments in **Comment** box.

Figure 6-159 Using LI-RADS



NOTE

The character limit for comment is 512 characters.

NOTE

Click the icon in the upper right corner of the LI-RADS(ACR) menu. A window will appear with reference information. Click the icon again to close the window.

Operator can also activate LI-RADS(ACR) by pressing **Report** on the Touch Panel.

1. Press **Report**, click **LI-RADS(ACR)** to activate LI-RADS(ACR).
2. Select **US category** and **US visualization score**, input comments in **Comment** box.

NOTE

There is no reference information when using Report button to activate LI-RADS.

LI-DARS(ACR) Report

After completion, press **Report** > **LI-RADS(ACR)** > **Whizz Report** to display the report, this applies to both system configurations.

US LI-RADS® Management

The table summarizes the suggested follow-up for patients based on the US LI-RADS® Category. Note that the US Visualization Score is not used to determine follow up.

Table 6-55 US LI-RADS Management

Category	Follow Up Action
US-1 Negative	Repeat surveillance US in 6 mo
US-2 Subthreshold	Repeat surveillance US in 3-6 mo ^{a,b}
US-3 Positive	Further characterization with multiphase contrast-enhanced imaging ^c
Footnotes: a. For observations < 10mm in size (US-2 Subthreshold), follow-up in 3-6 months is suggested. If observation does not show growth over a 2-year period, observation can be benign and study can be categorized as US-1 Negative. b. This is concordant with 2010 AASLD recommendations. c. Multiphase imaging may include multiphase contrast-enhanced CT, MR, or CEUS.	

Whizz Report

The ultrasound system enables the generation of patient reports based on the examination performed and the analyses that were made during the exam. The reports are generated using the data stored in the system with pre-selected template.

You may edit a report while performing the exam; customize, delete, or add measurements; and save changes until you use the Store command. Once Stored, the reports are read-only.

It is recommended that the data be saved often, and then carefully reviewed before the report is Stored. Use the worksheet to facilitate the review and adjust data before storing a report. The final report can be printed on a standard printer.

This ultrasound system allows users to create their own report template by a separate tool called **Whizz Report Editor** on users' personal computer.

The template can be saved into an USB and imported to the ultrasound system.

Creating a report

Reports summarize the data obtained in the examination. They can contain data, images, and cine loops.

The ultrasound system retrieves the patient information from the archive which is then used to generate the report. Once generated, the report can be viewed, images can be added, and the patient's personal data can be modified. The examination data itself CANNOT be changed.

Activating the Report

1. Select **Report > Whizz Report** on the Touch Panel.
2. Default report for the current application displays on the monitor.
The information entered during the examination is automatically filled in the appropriate fields (e.g. demographic, diagnosis, comments).
3. Use cursor to select images on the clipboard and add to the **Images** column.
4. Use report button controls on touch panel for more further action.

Table 6-56 Report Button Controls

Button	Description
Worksheet	Accesses Worksheet Page
Graph	Accesses OB Graph page (applies only to OB).
Anatomical Page	Accesses Anatomical Survey page (applies only to OB).
Template	Selects templates from the list of selected applications.
Store	Stores the report page into Archive as CHM file. NOTE: CHM is a compressed HTML help file.
Save As	Exports the report page to storage media as PDF format.

Button	Description
Print	Prints out the report to the default printer.
Retrieve	Retrieves the report page from Archive. Stored Date/Time is appended to the name of stored report.
Print Preview	Preview the report before printing.
Add Bull's Eye	Add Bull's Eye graph in report. (applies only to Cardiac).
Exit	Exit to scanning page.
Delete	Delete the report.

Selecting another template

You can select another template for the current patient:

1. Select **Report > Whizz Report > Template** on the Touch Panel.
2. A list of available templates and exam categories displays.
3. Select the desired template using the **Trackball** and press **Set** key.

The selected template displays on the monitor.

NOTE

If you select another exam category, the template list of the selected category displays. Select the desired template.

4. Select the desired template name and press Set.
5. The report changes to the selected template.

Factory Templates

The system has factory templates for each application. You can also create your own templates with **Whizz Report Editor** on a personal computer and import to the ultrasound system. You need to save revised/new templates with unique names.

A template may include one or more of the following:

- Measurements
- Worksheet or Vessel Summary Images
- Anatomical Surveys or Biophysical Profiles
- Anatomical Graphics
- Graphs
- Images areas
- Score Boxes

User-defined report templates can be added from the **Utility > Reports** menu. Refer to *Import Report Template to Ultrasound System* on page 580 for detail information.

Whizz Report Editor

Whizz Report Editor is a software application installed on user's personal computer. You can use it to design and create your own customized template from a blank template page, or you can use an existing template (factory or user-defined) and save the changes.


Installation and Operating environment

Whizz Report Editor only supports running on Win10 and above operating systems.

User can export the **Whizz Report Editor** from the ultrasound system to install on your own computer.

1. Press **Utility > Whizz Report** on the Touch Panel.
2. Click **Export Whizz Report Editor** icon on the monitor.

Figure 6-160 Export Whizz Report Editor



3. Select the removable media to save the **Whizz Report Editor** application.
4. Connect the removable media to your local computer with Win 10 and above operating system and double click the application to install.

Design or edit a report template

1. Double click the **Whizz Report Editor** application to open the tool on your personal computer.

2. The main page displays as below.

Figure 6-161 Whizz Report Editor main page




Table 6-57 Whizz Report Editor main page

No	Name	Description
a	Components	Lists all the components user may need for the report.
b	Tool Bar	List the tools for editing the main page of the template.
c	Components Editor	Use to edit the details of each selected component.
d	Template	Template preview for current design status.

Whizz Report Components

There are 4 categories of components for users to choose. User can choose general components under **General** and select special components from **OB**, **Carotid** and **Cardiac** for specific application.

Advanced Features

To add a component in the template, you can click and hold the cursor to drag the desired component to the report template.

Figure 6-162 Drag components to template




Table 6-58 Whizz Report Components

Category	Component	Description
General	Headers	This is the header of the report.
	Images	This is a box displays the images from the test.
	Conclusion	It's a input box to input conclusion in text format.
	Patient Information	This is a box displays patient information.
	Measurement	This is a box displays measurements.
	Referral Reason	It's a input box to input referral reason in text format.
	Diagnosis	It's a input box to input diagnosis in text format.
	Signature	It's a box to show Signature, user can write text in blank area.
	Anatomical Survey	This is a box to show anatomical survey result.
	Biophysical Profile	This is a box to show Biophysical Profile result.
	Line	Use to insert a line between two components.
	Measurements CV	This is a box displays measurements CV.
	Footers	This is the footer of report template, Operator and Perf. Physician can be changed to other patient information elements.
	Ratios Measurements	This is a box displays Ratios Measurements.
	Text Field	This is a Text input box, maximum number of this components is 20 in a report template.
	Fixed Text	This is a Fixed Text box, and user can edit text in designer tool, maximum number of this components is 20 in a report template.
OB	OB Graphs	This is a box displays OB Graphs. User can check or uncheck OB Graphs in worksheet->Graph.
	OB Calculations	This is a box to show calculations text. System will insert calculation measurements data to the table automatically.
	Growth Bar Graph	This is a box displays Growth Bar Graph. User can select fetus.
	OB Measurements	This is a box displays OB Measurements.
Carotid	Carotid Measurements	This is a box displays carotid measurements, In middle, this is a carotid image, and in left and right, there are some measurements, include measurement title and measurements data.

Category	Component	Description
Cardiac	Bull Eye Wall Motion Score Box	Displays two cardiac images and every area can be paint by different color.
	Cut Planes Wall	Displays two rows and per row has five cardiac images, and every area can be paint by different color.
	Score Table	This is a list to show Score and color of Score .

Tool Bar

On the top of the Whizz Report Editor main page, there is a tool bar.

Figure 6-163 Whizz Report Editor Tool Bar



Table 6-59 Whizz Report Editor Tool Bar

No	Tool Name	Description
1	Report Templates	Click to select report template from Default Templates or My Templates.
2	Report Display Ratio	Use drop-down menu to select report display ratio. Available ratio: 25%, 50%, 75%, 100%, 125% and 150%.
3	Undo/Redo	Click to Undo or Redo previous / next steps.
4	Move	Click on the icon, the icon will turn to blue and the cursor shape will change. You can move the template by click and hold the cursor. <i>NOTE: Press the icon again to disable the movement.</i>
5	Zoom Out/ Zoom In	Click to zoom out / zoom in the template.
6	Page Margin	Click to select template page margin from Normal, Narrow or Wide.
7	Template Name	Click the tab to input the template name. Default text is Default Template Name.
8	Export Html	Click to export current template in Html format.
9	Save as	Save current template as a new template.
10	Save	Save the template.

Whizz Report Component Configuration

On the right side of Whizz Report Editor main page, there is a component configuration column. You can configure the display for each component.

Private Options

Private Options is different for each component. Once user select a component in the template, all the configurable elements will be displayed under **Private Options**.


NOTE

Some components doesn't have configurable elements. Please skip if the Private Options drop-down menu is empty.

Headers Configuration

Header configuration allows user to edit report logo and hospital information.

Figure 6-164 Headers Configuration




1. **Logo** - Select desired logo for the report template.
2. **Delete** - Delete selected logo.
3. **Import** - Import new graphic for logo.
4. **Hospital Info** - Check the box and type text in the tab to show Hospital name, Address and Phone on the template.

Images Configuration

User can select the Grid to adjust the display of template images.


Figure 6-165 Images Configuration



Patient Information Configuration

This section allows user to configure patient information

Figure 6-166 Patient Information Configuration




1. **Show Title** - Enable or disable the component title.
2. **Table** - Select the design of the columns showing in the components.
3. **Info filter** - Check to select the information that displays under **Patient Information**.

Conclusion Configuration

Use can enable and disable the Summary in the component.

Figure 6-167 Conclusion Configuration




1. **Show Title** - Enable or disable the component title.
2. **Summary** - Display or hide **Summary** in the component.

Line Configuration

Use can change the **Stroke Width** of **Line** in the template.


Figure 6-168 Line Configuration



Measurement CV Configuration

Use can check to select desired **Vascular Measurements**.


Figure 6-169 Vascular Measurement Configuration



Footer Configuration

Use can check to select the information to display on footer.


Figure 6-170 Footer Configuration



Bull Eye Wall Motion Score Box Configuration

Use can select the Grid to adjust the display of template images.


Figure 6-171 Bull Eye Wall Motion Score Box Configuration



Container Options

Container Options allow users to customize the Padding, Stroke and Font for each component.

Figure 6-172 Container Options



1. **Padding** - Use to adjust the padding for the component.

2. **Stroke** - Use to adjust color, style and width for the component.
3. **Font** - Use to adjust color, size and style of the fonts in the component.

NOTE


Component title is not adjustable.

Report Template Preview

This section displays current editing template. User can preview during editing. You can move the components up and down or adjust the size of each component.


User can use the cursor to click and hold on the components to drag the component to move the component up and down.

Figure 6-173 Move the component up and down



User can move the cursor to the right bottom of each component till the cursor turns to arrow, click and hold to resize the component.

Figure 6-174 Resize the component



Export Report Template

User can export your own templated designed by **Whizz Report Editor** and import to ultrasound system.

Advanced Features

To export your own template from **Whizz Report Editor**:


- Click the **Export** icon in the tool bar of **Whizz Report Editor** to save the current editing template to a removable media.

Figure 6-175 Export current template as Html Document



- Click **Report Templates** icon in the tool bar of **Whizz Report Editor**, then select and save desired template from **My Templates** > **Export**.

Figure 6-176 Export Html Documents from My Templates



Import Report Template to Ultrasound System

User can import your own templated designed by **Whizz Report Editor** to ultrasound system.

To import a new template from **Whizz Report Editor**:

- Save and export your own template to a removable media. Refer to *Export Report Template* on page 579 for exporting details.
- Connect the removable media to ultrasound system.
- Press **Utility** > **Whizz Report** on Touch Panel.
- Click **Import** to open the removable media.
- Select the template from the removable media and click **Import** to complete.

Once the import is completed, a templates successfully imported information window will pop out.


Click **OK** to continue.

6. Select the desired template from **Available Templates > My templates**, then click



to move the template under **Section**.

Figure 6-177 Import report template



a: **Set Default** - Once the imported template is selected, there Set Default icon displays. User can click to set imported template as default.

b: **Up and Down icon** - Use these buttons to move the selected template up and down.

c: **Reset**: Reset to original settings.

7. Press **Report > Whizz Report > Template** on the Touch Panel to change the template by **Trackball** and **Cursor** key.

Figure 6-178 Change the report template



Scan Coach

Scan Coach Manager

Before performing Scan Coach function on the system, please read and accept the Statement, Disclaimer and Limitation of liability described as below:

Statement

- Scan Coach IS NOT MEANT TO REPLACE TRAINING, OR TUTORIALS/HANDS-ON. IT IS MERELY A REFRESHER TOOL OF ALREADY RECEIVED EDUCATION AND TRAINING.
- Scan Coach IS AN ON -DEMAND REFRESHER AND REMINDER TOOL WHICH DISPLAYS INFORMATION IN THE FORM OF DICOM IMAGES, ANIMATIONS AND VOICE COMMENT; THE PROVIDED REFERENCE MATERIAL MAY HELP USER IN ACQUIRING ULTRASOUND IMAGES.
- Scan Coach DOES NOT PROVIDE A DIAGNOSIS, BUT DOES PROVIDE REFERENCE MATERIAL FOR SOME TYPES OF ACQUISITIONS.
- Scan Coach IS MEANT TO PROVIDE REFERENCE MATERIAL FOR ACQUISITION, BUT IT IS NOT INTENDED TO IDENTIFY DIAGNOSTIC IMAGE QUALITY. ACTUAL IMAGES, INCLUDING IMAGE QUALITY, OF THE DEVICE MAY VARY VERSUS THE PROVIDED REFERENCE MATERIAL.

Disclaimer

Scan Coach AND THE CONTENT AVAILABLE THROUGH IT ARE PROVIDED ON AN "AS IS" AND "AS AVAILABLE" BASIS. YOU EXPRESSLY AGREE THAT USE OF Scan Coach AND/OR ITS CONTENT IS AT YOUR SOLE RISK. TO THE FULLEST EXTENT PERMISSIBLE PURSUANT TO APPLICABLE LAW, GE AND ITS AFFILIATES DISCLAIM ALL WARRANTIES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR NON-INFRINGEMENT. YOU EXPRESSLY AGREE THAT USE OF Scan Coach, INCLUDING ALL CONTENT, IS AT YOUR SOLE RISK.

ANY RESPONSIBILITY FOR IMAGE ACQUISITION, IMAGE INTERPRETATION, IDENTIFICATION OF ANATOMICAL PARTS, ANATOMICAL MEASUREMENTS AND CLINICAL DIAGNOSIS LIES WITH YOU, THE ULTRASOUND USER.

TO THE FULLEST EXTENT PERMISSIBLE PURSUANT TO APPLICABLE LAW, GE AND ITS AFFILIATES DISCLAIM ALL LIABILITY AND RESPONSIBILITY FOR ANY DAMAGES SUFFERED BY A USER, OR PATIENT OF A USER INCLUDING, BUT NOT LIMITED TO, DAMAGES ARISING OUT OF LOSS OR INACCURACY OF CONTENT OF Scan Coach.

IT IS THE RESPONSIBILITY OF YOU, THE ULTRASOUND USER, TO UNDERTAKE ALL NECESSARY AND CUSTOMARY TRAINING, AND TO MEET ALL OTHER REQUIREMENTS OF APPLICABLE LAWS AND REGULATIONS, BEFORE USING THE ULTRASOUND DEVICE.

Limitation of Liability

UNDER NO CIRCUMSTANCES, INCLUDING, WITHOUT LIMITATION, NEGLIGENCE, SHALL GE OR ITS PARENTS, SUBSIDIARIES, AFFILIATES, OFFICERS, DIRECTORS,

EMPLOYEES, AGENTS, OR SUPPLIERS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES ARISING FROM OR IN CONNECTION WITH THE USE OF OR THE INABILITY TO USE Scan Coach.


Introduction

Scan Coach is based on Scan Assistant, it displays information which helps user acquire the right scan plane. The reference image indicates how the scan plane image for a given step should look like. The probe/beam anime image shows the corresponding probe placement, or beam formation for getting the correct scan plane. The schema anime shows the key anatomical structures to be visualized in two dimensional mode.

Scan Coach can be used as an on-demand refresher/reminder tool in live scanning.

Scan Coach Description

Figure 6-179 Scan Coach Display Description



1. Program step number and step name.

NOTE

The user is also able to select the program step by pressing the Up/Down key on the AN keyboard.

2. **List View** - Select this icon to display program list.
3. **Image View** - Select this icon to activate Animation and Reference image.
4. **Both View** - Select this icon to display both view of Scan Coach.
5. Select this icon to activate Scan Coach Edit page.

Advanced Features

6. Animation of probe position.
7. Two dimensional schema of the current scan plane.
8. Voice Comment list.
9. Recorder for voice comment.


Setting up Scan Coach

Go to **Utility > Scan Coach > Scan Coach Manager** to create, import/export and manage the Scan Coach programs.

Using Scan Coach with Voice Comment

1. Go to **Utility > Scan Coach > Scan Coach Manager**, select the protocol from the left column, move it to the **Program Selections** column. Then select **Exit** to the Scanning page.

Figure 6-180 Select Scan Coach Protocol



2. Press **Assistant** key on the control panel and then select **Scan Coach** on the Touch Panel to activate the Scan Coach function. The Scan Coach addendum displays. Select **Acknowledge** to continue.
Or select **Cancel** not to start Scan Coach.
Check the box before **Don't show next time** and the Scan Coach addendum will not display next time.
3. Select the protocol from the pull-down menu, and then select **Start**.


NOTE

If there is only one program in the current application, the system will display the Scan Coach screen directly. This step is only available when there are two or more programs.

4. Select program, for example, LIVER LENGTH, then select **Image View** icon to review the reference image and the probe position and schema.

It displays reference image, probe position, voice comment and schema to guide the user acquire the right scan plane.

Figure 6-181 Using Scan Coach 2



- a. live scan image
 - b. reference image
 - c. probe position
 - d. two dimensional schema of the current scan plane
 - e. voice comment
5. You can also press the **List View** icon to return to the program step list and select another program.

Using Voice Comment


Voice Comment is a quick method to record comments while Scan Coach.

To use the voice comment:

1. Connect the microphone to ultrasound system.

- Record the voice comments and save to **Voice Comment** list.

Figure 6-182 Voice Comment Recorder



- Click to start recording.
- Click to play the recorded comment.
- Click to stop recording.

NOTE

The maximum recording time is 2 minutes. Once the recording time reaches 2 minutes, the recording will be stopped automatically.

- Input the file name.
- Click to save the Voice Comment.


- All the recorded voice comments will be saved to the **Voice Comment** list.

NOTE

User can save maximum 5 voice comments, for more comments, please delete the existing comments to release space.

- Play the voice comments.

Figure 6-183 Voice Comment List



- a. Click to switch from playing each comments continuously on the list or playing selected comments only.
- b. Click to delete the voice comment.
- c. Click to play the voice comment.


Edit Scan Coach Protocol

1. To edit Scan Coach protocols, go to **Utility > Scan Coach > Scan Coach Manager** and select **Creator**.

Or Select the icon  in the program list. Then select  and the system will lead to the Scan Coach Edit page.


2. The Scan Coach Creator is used to build customized programs that can be imported onto the Versana Premier/Versana Premier Lotus. Select **Scan Coach** on the tool bar to display the Scan Coach edit page.

Figure 6-184 Scan Coach Creator




3. Select the Scan Plane. Select **Browse** to upload the reference image.

Figure 6-185 Select Scan Plane



4. Select **Upload** to upload the reference image from the folder called "Reference Images" of the corresponding application. Then select **Close** after the upload is complete.

Figure 6-186 Upload Reference Image



NOTE

The reference image only supports *.dcm file.

The reference image can also be uploaded from external USB stick/HDD/CD/DVD. Select the correct directory from the pull-down menu, and then select the appropriate image to upload.

5. Upload Probe Position Image and Schema Image from the folder called "Other Images" of the corresponding application, or from the external devices.

NOTE

A protocol can be created without uploading the reference image, but the protocol without reference image do not display in **Utility > Scan Coach > Scan Coach Manager**. It can be found via **Utility > Scan Coach > Scan Coach Manager > Custom Programs**.


Activating the Scan Coach Creator

To activate the Scan Coach Creator on the Versana Premier/Versana Premier Lotus, press **Utility > Scan Coach > Scan Coach Manager** and select **Creator**.

Scan Coach Creator

The Scan Coach Creator is used to build customized Programs that can be imported onto this ultrasound system. These Programs automate many of the steps normally performed manually by the user, thereby reducing the number of user actions and the amount of time to perform an exam.

Figure 6-187 Tool Layout Overview



1. Menu and Tool Bars
2. Scan Coach Button
3. Steps
4. Scan Coach Attributes
5. Rule Checking

NOTE

Scan Coach Button is not available for Scan Assistant, only for Scan Coach.

The Scan Coach Creator tool can be used both on the scanner and as an off scanner tool. Where there are differences in behavior, this user's guide uses the term 'on scanner' to indicate when the tool is running on the scanner and 'off scanner' to indicate when the tool is running off the scanner.

File Handling

When using Scan Coach Creator off the scanner, it is very important to organize the programs in a way that will make it easy to import the programs onto the scanner. Each Program is a computer file. While these computer files can be copied, pasted and deleted like any other computer file, the Program files are only viewable using the Scan Coach Creator.

Off-Scanner Directory Structure

The Scan Coach Creator organizes the Programs in a directory structure that allows easy importing into the Versana Premier/Versana Premier Lotus. In order to be imported, all Programs must be stored in a LOGIQ_SCAN_COACH Programs Directory. Within this directory, one or more user-specified directories are created. Within each of these user-specified directories are the category directories (VAS, ABD, etc.) that hold the actual Programs.

The dialog in the figure below allows the user to specify the location of the LOGIQ_SCAN_COACH directory (root directory) and to either select an existing User Program Directory or create a new one.

The Directory Structure dialog can be accessed via the File menu and File Toolbar.

File Extensions

Factory defined Programs have an .ep (exam Program) extension while user-defined Programs have an .uep (user exam Program) extension. Both factory and user-defined Programs can be read into the Scan Coach Creator, but only user-defined Programs are created. If a factory Program is read into the Scan Coach Creator and then edited, it is saved as a user-defined Program.

Upon installation of the Scan Coach Creator, files with a .ep or .uep extension are automatically associated with the Scan Coach Creator.

Exporting Programs from Ultrasound System

Factory or user-defined Programs on the ultrasound system are easily exported for editing with the Scan Coach Creator.

On the ultrasound system:

1. Insert a USB storage device (or CD/DVD).
2. Select **Utility > Scan Coach**.
3. Select **Export**.
4. Select the media type and specify a directory. If a directory is specified that already exists, the Export adds the Programs along with any existing Programs. If the names of Programs are the same, use the resulting dialog to decide how to continue.
5. Select the Program to be exported and export them.

Importing Programs to Ultrasound System

Programs created with the Scan Coach Creator are easily imported to the Versana Premier/Versana Premier Lotus.

Copy the complete LOGIQ_SCAN_COACH directory from the system to a USB device (or CD/DVD). The LOGIQ_SCAN_COACH directory needs to be at the top level (not in a subdirectory) on the USB device (or CD/DVD).

Example of Directory structure:

LOGIQ_SCAN_COACH

MyUserNameDirectory

ABD

CARD

GYN

OB

VAS

On the Versana Premier/Versana Premier Lotus:

1. Insert the USB device (or CD/DVD).
2. Select **Utility > Scan Coach**.
3. Select **Import**.
4. Select the media type.
5. Select the Programs to be imported and import them. If you attempt to import Programs that already exist with the same name, use the resulting dialog to decide how to continue.

Creating New Programs

A new Program is created by selecting **File > New**, by clicking on the New document icon in the Toolbar, or by using the keyboard shortcut Ctrl+N.

1. Before creating a New Program, select **Single Step** in the Toolbar.

Figure 6-188 File Toolbar




2. Proceed to add/update your settings for the Step: Step Name, Instructions, etc.
3. Once finished, highlight the finished Step.

Figure 6-189 Highlight Step to Copy



4. Select **Edit > Copy**.

Figure 6-190 Edit -> Copy



5. In the Toolbar along the left, select **Insert Step Before Selected** or **Insert Step After Selected**.

Figure 6-191 Insert Step



6. Highlight the copied step and proceed to edit accordingly.
7. Proceed to follow the same procedure to add more steps to your Program.
8. When you are done, select **Check** to verify your Steps.
9. The results are listed as to whether the Scan Coach Rule Check Passed or if any Issues were detected. Issues found when running the check do not mean the Program is unusable.

NOTE

The rule check may report an unequal number of left and right steps. This may or may not be the expected result. If a change is made in response to the rule check results, a new rule check can be run to see if the issue has been resolved.

Views

A Program is made up of a series of steps. Each step is made up of various step attributes. The step and step attribute data can be viewed in many ways using the Scan Coach Creator.

The different ways to look at the data are called Views. The view of choice is selected from the **View** Selection Menu or the View Toolbar Menu.

Figure 6-192 View Selection Menu and View Toolbar Menu



Single Step Views

There are two Single Step views: Basic and All. The Basic view shows the most common attributes of the selected step. The All view shows all of the attributes of a given step.


For both views, the step names are shown on the left with the active step highlighted. The step attributes appear on the right and are separated into four groupings:

General attributes at the top

Imaging and Comment attributes on the left

Measure attributes on the right.

Figure 6-193 Basic Single Step View



1. Steps
2. Step Attributes

Scan Coach Features

Scan Coach allows the user to program the steps in an exam and to program certain attributes for each step. The attributes are what give the Scan Coach Program behavior. The tables below provide the names of all attributes along with a description of how each one is interpreted by the Scan Coach feature.

General Attributes

Table 6-60 General Attributes

Attribute Name	Selections	Description	First Step Default	Other Steps Default
Store Order	(Not Applicable)	Button used to enable the Store Order Definition dialog so that the Store Order can be set		
Advance On	Store	Advance to the next step and go live after Print / Image Store (e.g. P1 key). This can be a single image store or a loop store.	Print	Same as previous step
	Store & Unfreeze	Advance to the next step after Print / Image Store (e.g. P1 key) and unfreeze. This can be a single image store or a loop store.		
	User Selection	Advance to next step only after next step is manually selected (e.g. down arrow)		
Instructions	[Any Text]	User notes displayed in the Scan Coach Navigation menu when the step is active	Blank	Same as previous step
Optional	Optional (checked)	An optional step is given a check mark during Program execution even if no image is acquired	Mandatory	Same as previous step
	Mandatory (unchecked)	A mandatory step is given a check mark only if an image is acquired for the step		
Learn Probe	On (checked)	Learn and change the probe for the user	Off	Off
	Off (unchecked)	No probe change		

Comment Attributes**Table 6-61 Comment Attributes**

Attribute Name	Selections	Description	First Step Default	Other Steps Default
Comment 1, Comment 2	[Any Text]	User annotation associated with the step. When editing in a Multi Step View, use Alt+Enter to create a new line.	Same as Step Name attribute	Same as Step Name attribute

Attribute Name	Selections	Description	First Step Default	Other Steps Default
Location 1, Location 2	Top Left	Annotation is placed in the top left corner of the image area	Location 1: Bottom Center	Same as previous step
	Middle Left	Annotation is placed in the middle left side of the image area		
	Bottom Left	Annotation is placed in the bottom left corner of the image area	Location 2: Dual Right: Bottom Center	
	Top Center	Annotation is placed in the top center of the image area		
	Bottom Center	Annotation is placed in the bottom center of the image area		
	Top Right	Annotation is placed in the top right corner of the image area		
	Mid Right	Annotation is placed in the middle right side of the image area		
	Bottom Right	Annotation is placed in the bottom right corner of the image area		
	Dual Left: Bottom Center	Annotation is placed on the bottom center of the left image in dual screen		
	Dual Left: Top Left	Annotation is placed on the top left of the left image in dual screen		
	Dual Left: Bottom Right	Annotation is placed on the bottom right of the left image in dual screen		
	Dual Right: Bottom Center	Annotation is placed on the bottom center of the right image in dual screen		
	Dual Right: Top Left	Annotation is placed on the top left of the right image in dual screen		
	Dual Right: Bottom Right	Annotation is placed on the bottom right of the right image in dual screen		

Attribute Name	Selections	Description	First Step Default	Other Steps Default
BP1, 2	Blank	Body Pattern not specified. Scan Coach does not set Body Pattern.	Blank	Same as previous step
	Body pattern graphics with or without probe position	Selected Body Pattern with or without probe position will be set.	Blank	Same as previous step
BP Specify	(Not applicable)	Button used to enable the Body Pattern Selection dialog so that the Body Pattern graphic can be selected and probe position can be set		
BP Clear	(Not applicable)	Clears BP1, BP2 defined for the step		
BP Probe	On (checked)	BP Probe mark will be set by Scan Coach	Off	Same as previous step
	Off (unchecked)	Probe mark not specified. Scan Coach does not set Probe mark.		

Imaging Mode Attributes

The probe and application associated with a program is not configurable. Instead, the scanner remembers the last probe and application used for a given Scan Coach program and automatically selects them the next time the program is started.

Table 6-62 Imaging Mode Attributes

Attribute Name	Selections	Description	First Step Default	Other Steps Default
Color Mode	On (checked)	Color Doppler is on	Off	Same as previous step
	Off (unchecked)	Color Doppler is off		
PW Doppler Mode	On (checked)	PW Doppler is on. If PW Doppler is not on and the new activated step indicates that Doppler should be on, the Mode Cursor is displayed or, if the Mode Cursor is already displayed, then PW Doppler is turned on.	Off	Same as previous step
	Off (unchecked)	PW Doppler is off		

Attribute Name	Selections	Description	First Step Default	Other Steps Default
M Mode	On (checked)	M-Mode is on	Off	Same as previous step
	Off (unchecked)	M-mode is off		
CW Mode	On (checked)	CW Mode is on	Off	Same as previous step
	Off (unchecked)	CW Mode is off		
Dual	Off	Dual screen not in use	Off	Same as previous step
	Left Active	Dual screen is active and the left image is the active image.		
	Right Active	Dual screen is active and the right image is the active image.		
	DualView (simul)	DualView is active (both left and right images are live)		
	Dual on Freeze	DualView is active on freeze		
3D	On (checked)	3D on	Off	Same as previous step
	Off (unchecked)	3D off		
LOGIQView	On (checked)	LOGIQView on	Off	Same as previous step
	Off (unchecked)	LOGIQView off		

Imaging Preference Attributes

Imaging Preferences work slightly different than other attributes. For example, if an abdomen Program has 20 steps and all steps have the Harmonics attribute set to Default, then Scan Coach will not affect the harmonics setting. Now, assume that steps 10-12 are gallbladder steps and that the harmonics attribute has been set to on for these steps. When transitioning into this group of steps (step 9 to step 10, e.g.), harmonics will be turned on (or remain on if it was previously on). If harmonics is then manually turned off in step 10 then Scan Coach will not turn it back on when advancing to step 11. In other words, a group of consecutive steps with the same Imaging Preference are treated as a group by Scan Coach and not as individual steps.

Table 6-63 Imaging Preference Attributes

Attribute Name	Selections	Description	First Step Default	Other Steps Default
Harmonics	Off	Harmonics off	Default	Same as previous step
	On	Harmonics on		
	Default	Harmonics not specified. Scan Coach does not set Harmonics on or off.		
Virtual Convex	Off	Virtual Convex off	Default	Same as previous step
	On	Virtual Convex on		
	Default	Virtual Convex not specified. Scan Coach does not set Virtual Convex on or off.		
CrossXBeam	Off	CrossXBeam off	Default	Same as previous step
	On	CrossXBeam on		
	Default	CrossXBeam not specified. Scan Coach does not set CrossXBeam on or off.		
PDI	Off	PDI off	Default	Same as previous step
	On	PDI on		
	Default	PDI not specified. Scan Coach does not set PDI on or off.		
Color/Doppler Steer	Left	Color/Doppler steered to the left	Center	Same as previous step
	Center	Color/Doppler not steered		
	Right	Color/Doppler steered to the right		

Measurement Attributes
Table 6-64 Measure Attributes

Attribute Name	Selections	Description	First Step Default	Other Steps Default
B/M Measure Trigger	Measure Key	Initiate "Measure 1" attribute when the Measure key is manually selected.	None	None
	Freeze Key	Initiate "Measure 1" attribute when the image is frozen.		
	Store Image	Initiate "Measure 1" attribute when the Measure key is manually selected or the image is stored. This is used to store / print an image and then measure on it and then store it again. Therefore, the Advance On Print attribute is ignored on the first store / print when the Measure Trigger attribute is set to Image Store.		
	None	Measurements are not triggered by Scan Coach. The "Measure 1" attribute is ignored.		
Side	Rt	The side measurement qualifier is set to Right side of the body	Rt	Derived from Step Name attribute if possible. Otherwise, same as previous step
	Lt	The side measurement qualifier is set to Right side of the body		
	None	The side measurement qualifier is not used (neither Right nor Left)		

Attribute Name	Selections	Description	First Step Default	Other Steps Default
Fetus	A	The fetus measurement qualifier is set to Fetus A	A	Same as previous step
	B	The fetus measurement qualifier is set to Fetus B		
	C	The fetus measurement qualifier is set to Fetus C		
	D	The fetus measurement qualifier is set to Fetus D		
Location	Prox	The location measurement qualifier is set to Proximal	Prox	Same as previous step
	Mid	The location measurement qualifier is set to Middle		
	Dist	The location measurement qualifier is set to Distal		
	None	The location measurement qualifier is not used		
Measure 1	Various 2D or M-mode measurements	Specifies the first 2D or M-mode measurement to be initiated. The point at which the measurement is initiated is based upon the Measure Trigger attribute. See Section 6 for more information.	Blank	Blank
Measure 2	Various 2D or M-mode measurements	Specifies the second 2D or M-mode measurement to be initiated after the measurement associated with the Measure 1 attribute is completed. See Section 6 for more information	Blank	Blank
Measure 3	Various 2D or M-mode measurements	Specifies the second 2D or M-mode measurement to be initiated after the measurement associated with the Measure 2 attribute is completed. See Section 6 for more information.	Blank	Blank
Dop Vessel Assign	Various Doppler measurement Vessel folders	Specifies the Vessel folder to assign auto calcs to. The assignment happens when the image is stored / printed (P1 key, e.g.). See Section 6 for more information.	Blank	Blank

Attribute Name	Selections	Description	First Step Default	Other Steps Default
Auto Calcs	Default	Auto Calcs state not specified. Scan Coach does not set Auto Calcs state.	Default	Same as previous step
	Off	Auto Calcs state set to Off		
	Frozen	Auto Calcs state set to Frozen		
	Live	Auto Calcs state set to Live		
Auto Calc Params	Various Auto Calc parameters	Specifies the auto calc parameters to be used.	Default	Same as previous step
	Default	Auto Calc parameters are not specified. Scan Coach does not set the Auto Calc parameters.		
Auto Calc Specify	[Not Applicable]	Button used to enable the Auto Calcs Parameter Selection dialog so that the Auto Calc Params attribute can be set	[Not Applicable]	[Not Applicable]
Auto Calc Default	[Not Applicable]	Button used to set the Auto Calcs Params attribute to Default.	[Not Applicable]	[Not Applicable]
Double Print	On (checked)	If an Image Store / Print (P1 key, e.g.) is performed on an image with measurements, the image is stored / printed two times, once with the measurements and once without.	Off	Same as previous step
	Off (unchecked)	No special Store / Print behavior.		

Attribute Name	Selections	Description	First Step Default	Other Steps Default
Double Print	On (checked)	If an Image Store / Print (P1 key, e.g.) is performed on an image with color, the image is stored / printed two times, once with color and once without. If double print on color and double print on measurements are both configured to be on, the image is stored / printed two times, once with the measurements and once without.		
	Off (unchecked)	No special Store / Print behavior.		
View Worksheet	On (checked)	The worksheet is turned on	Off	Off
	Off (unchecked)	The worksheet is not turned on		

Numerics

- 3D Acquisition 396
- 3D/4D imaging
 - introduction 508
 - manipulating the volume of interest 525
 - operational controls 515
 - principles of operation 508

A

- accessories
 - ordering 250
 - requesting a catalog 250
- ALARA (as low as reasonably achievable), bioeffects 260
- Allergic reactions to latex-containing medical devices 269
- Alpha Hip, pediatric measurement 142
- area measurements
 - ellipse 97
 - spline 98
 - trace 97, 98
- Auto Vascular Calculation, see also Manual Vascular calculations
 - activating 131
 - setting up calculation parameters 130
- Auto vs. Manual Calculations 143

B

- B-Mode measurements
 - Echo Level 100
- B-Mode measurements, general 96
- B-Mode measurements, generic
 - Volume 104
- B-Mode measurements, mode
 - circumference and area (ellipse) 97
 - circumference and area (spline trace) 98
 - circumference and area (trace) 98
 - distance 96
- biological hazards 264
- Bladder Volume 135
- BluetoothSave 168

C

- calculations
 - urology 134
- Caution icon, defined 260
- Check System Date and Time 46
- CINE loop
 - previewing only 163
 - storing and previewing 162
 - storing without previewing 163
- circumference measurements
 - ellipse 97
 - spline 98
 - trace 97, 98
- Connectivity Functions 186
- contacts
 - clinical questions 250
 - Internet 250
 - service questions 250
- contraindications 14
- Control panel adjustment 29
- Curved Anatomical M-Mode 82

D

- Danger icon, defined 260
- Database Protection Mechanism 48
- dD Ratio, pediatric measurement 143
- Disposal 205
- distance measurement
 - general 96
- Doppler measurements, mode
 - TAMAX and TAMEAN 102
 - time interval 104
 - velocity 102
- Doppler Mode, general measurements 102

E

- E7C8L-RS Operation 316
- ECG Cleaning and Disinfection 386
- ECG Operation 385
- ECG Patient preparation 386
- Echo level measurements 100
- Elastography 404

Index

Elastography, 2D measurement 406
electrical hazard 264
electromagnetic compatibility (EMC) 270
ellipse measurement, general 97
EMC (electromagnetic compatibility) 270
Endometrium thickness (Endo), GYN exam 118
equipment safety 264
exam
 OB 109
explosion hazard 264
Export Report Template 579

F

Federal law (USA), requirements 15
Fetal growth curve graph
 multiple fetus 115
Fetal trending
 multiple fetuses 115
Fetus
 entering number of 114
Fetus Compare
 multiple fetus 115
Follicle measurements, GYN exam 118

G

Gel Warmer (Option) 26
General imaging
 changing presets 155
GYN exam
 endometrium thickness 118
 follicle measurements 118
 ovaries measurements 118
 uterus measurements 119
Gynecology exam 117

H

hazards, safety symbols 260
hazards, types
 biological 264
 electrical 263, 264
 explosion 264
 mechanical 263
 smoke and fire 264
Hip Dysplasia, pediatric measurement 141

I

Imaging Insights 187
Imaging presets, changing
 General 155
Import Report Template 580
Indicator Light 30
information, requesting 250

L

Loading the Software 222
Log Off 47

M

M-Mode
 CAMM 82
M-Mode measurements, mode
 time interval 101
 tissue depth 101
M-Mode, general measurements 100
Manual Vascular Calculation 133
Multiple fetuses 114
 fetal trending 115
 identifying 115
 on OB worksheet 116
My Trainer 468

N

Noisy ECG Signal 388

O

OB
 graph 113
 identifying multiple fetuses 115
OB exam
 preparing 109
OB graphs
 patient data 114
 viewing 113
OB multigestational 114
OB worksheet
 multiple fetuses 116
Optimizing M Color Flow 86
Orientation Help, 3D/4D imaging 518
ovaries, measurement 118

P

- Patient data
 - OB graphs 114
- patient safety 262
- pediatric measurements, types
 - alpha HIP 142
 - dD ratio 143
 - hip dysplasia 141
- Phantoms 234
- Physical A/N keyboard 30
- prescription device, caution 15
- Preset Program menu
 - Acoustic Output
 - Fetal Exposure 109
 - Prudent Use 109
- probe orientation, 3D/4D imaging 514
- Prostate Volume 136
- prudent use 260

Q

- QAnalysis 476
 - Drift Compensation 484
 - Exit 487
 - generating a trace 479
 - manipulating the sample area 480
 - Select Image Range 478
 - smoothing 484
- QAPlot Control
 - horizontal sweep 484
- Quality Assurance
 - baselines 234
 - frequency of tests 233
 - periodic checks 235
 - phantoms 234
 - record keeping 241
 - system setup 235
 - test descriptions 236
 - typical tests 233

R

- Record keeping 241
- reference images, 3D/4D 512
- Renal Volume 135
- Render view, 3D/4D imaging 517
- Report Templates 569
- Reporting 262

S

- safety
 - electromagnetic compatibility (EMC) 270
 - equipment 264
 - hazards 260, 264
 - patient 262
 - electrical hazards 263
 - mechanical hazards 263
 - patient identification 262
 - personnel 264
 - precaution icons, defined 260
 - precaution levels, defined 260
- Scalpel, 3D/4D imaging 528
- Scan Assistant
 - Availability 460
 - Definitions 460
 - Description 461
 - Setting Up 461
 - Using 464
- Scan Coach Creator 588
- Scan Coach Manager 582
- Scan Coach Protocol 587
- Scan Coach with Voice Comment 585
- Sectional view, 3D/4D imaging 517
- service, requesting 250
- Setting up Scan Coach 584
- Small Parts 108
- Static 3D imaging 526

T

- TAMAX and TAMEAN, Doppler mode
 - measurement
 - manual trace 102
- TI and MI Display Accuracy 283
- TIC analysis
 - activating 414
 - curve fitting 430
 - exiting 416
 - generating a trace 419
 - manipulating the sample area 422
 - smoothing 427
- TIC Plot Control
 - horizontal sweep 427
 - line style 427
 - overview 425
 - vertical auto scaling 426
- Time interval
 - Doppler mode measurement 104

Index

Time interval ()

M-Mode measurement 101

Tissue depth, M-Mode measurement 101

Trace measurement, general 97, 98

TVI

Quantitative Analysis 476

U

urology measurements

pelvic floor 138

Using Whizz Follicle 537

Uterus measurements 119

V

Velocity, Doppler measurement 102

Volume measurement, B-Mode generic
measurement 104

W

Warning icon, defined 260

Whizz Report 568

Whizz Report Editor 570

Whizz Report Tool Bar 574

