ANIOSPRAY QUICK

Disinfectant
for spraying

- Ready to use solution
- Ethanol-based composition
- Broad activity spectrum in short contact time
- Does not leave any trace when drying
- Pleasantly fragranced
- Aldehyde free, without colouring

INDICATIONS

Rapid action disinfection of previously cleaned, non-immersible medical devices resistant to alcohol (stethoscopes, pressure sensors, blood sugar testers, etc.).

CHARACTERISTICS

- Ready to use product
- Clear colourless solution
- Density : approx 0.9
- pH at +20°C : approx 5.5

30 s.
mini


READY
TO USE

ANIOSPRAY QUICK

Disinfectant for spraying

INSTRUCTIONS FOR USE


To be used with manual sprayer ANIOS equipment.

1 Previously clean the medical devices with an ANIOS detergent disinfectant.

2 Spray an even film on the medical devices surfaces.

3 Leave for a **30 second minimum** contact time.

COMPOSITION

Hydroalcoholic solution (ethanol 55%), Quaternary ammonium propionate, perfume.

PRECAUTIONS FOR USE

Dangerous – respect the precautions for use (Drawn up according to the European rules in force regarding the classification and labelling of chemical products).

Store between +5°C and +25°C.

Class IIa medical device (Directive 93/42/EEC as amended).

PACKAGING

1 12 bottles, 1L each with 4 sprayers Ref. 2084.073

2 4 cans, 5 litres each Ref. 2084.034

MICROBIOLOGICAL PROPERTIES

Active against	Standards	Contact time
Bacteria	EN 1040, EN 13727, EN 14561, EN 13697. SARM (EN 14561)	30 sec.
Mycobacteria	Mycobacterium terrae : EN 14348, EN 14563	30 sec.
	Mycobacterium avium : EN 14348, EN 14563	5 min.
Yeasts / Moulds	Candida albicans : EN 1275, EN 13624, EN 14562, EN 13697. Tricophyton mentagrophytes : EN 14562	30 sec.
	Aspergillus niger : EN 1275, EN 13624, EN 14562, EN 13697.	
Viruses	Adenovirus, HIV-1, PRV (surrogate of HBV), BVDV (surrogate of HCV), Vaccinia virus* Herpes-virus, Coronavirus, Rotavirus, Influenza virus A (H ₁ N ₁)	30 sec.

* RKI : The Robert Koch Institute recommends the use of disinfectants claiming "begrenzt viruzid wirksan", i.e. limited virucidal activity to interrupt the infection-chain. According to the RKI, the limited virucidal activity corresponds to an activity against vacciniavirus and BVDV.

All the product's antimicrobial activity is included in the scientific dossier, available on request

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IMDG : International Maritime Dangerous Goods.

IATA : International Air Transport Association.

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RID : Regulations concerning the International carriage of Dangerous goods by rail.

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DECLARATION CE DE CONFORMITE DECLARATION OF CONFORMITY	Revision: 1
	Page: 1 of 1

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respecte toutes les dispositions de la directive 93/42/CEE.
meets all the provisions of the Directive 93/42/EEC.

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ANNEX II excluding section 4 Directive 93/42/EEC concerning medical devices

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For class III devices, a EC design certificate is required

Fabricant / Manufacturer / Fabricante

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Detergent disinfectant or disinfectant for invasive and/or non invasive devices

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GMED certifies that, on the basis of the results contained in the file referenced P177315, the quality system - for design, manufacturing, and final inspection - of medical devices listed here above complies with the requirements of the Directive 93/42/EEC, annex II excluding section 4.

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The validity of the certificate is subject to periodic or unexpected verification

Début de validité / Effective date / Fecha efectiva : November 20th, 2019 (included)

Valable jusqu'au / Expiry date / Fecha de expiración : June 10th, 2021 (included)

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Detergents disinfectants or disinfectants for invasive and/or non-invasive devices :

- Disinfectants for dialysis circuits
- Disinfectants for manual use and/or automatic machines for invasive devices
- Detergent disinfectant for manual use and/or automatic machines for invasive devices
- Detergent disinfectant for manual use and/or automatic machines for non-invasive devices
- Detergent disinfectant or disinfectant sprays for invasive devices
- Detergent disinfectant or disinfectant sprays for non-invasive devices
- Detergent disinfectant or disinfectant wipes for invasive devices
- Detergent disinfectant or disinfectant wipes for non-invasive devices

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GMED 0459

**On behalf of the President
Béatrice LYS
Technical Director**

Identification des dispositifs / Identification of devices / Identificación de los dispositivos

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Medical devices covered by this certificate are referenced on the manufacturer's list of products dated november 19th 2019 (13 pages) authenticated by GMED on november 20th, 2019

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GMED 0459

On behalf of the President
Béatrice LYS
Technical Director

ATTESTATION CE / EC CERTIFICATE / CERTIFICADO CE

Approbation du Système Complet d'assurance Qualité / Approval of full Quality Assurance System / Aprobación del sistema completo de Seguro de la calidad

ANNEXE II excluant le point 4 Directive 93/42/CEE relative aux dispositifs médicaux / Annex II excluding section 4 directive 93/42/EEC concerning medical devices / ANEXO II excluyendo el punto 4 directiva 93/42/CEE relativa a los productos sanitarios
Pour les dispositifs de classe III, un certificat CE de conception est requis / For class III devices, a EC design certificate is required.

Fabricant / Manufacturer / Fabricante

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*Design, development, manufacturing, servicing and sales of antimicrobial products, detergents and materials
for disinfecting processes of medical devices dedicated to hospitals and clinical sectors.*

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