



TRASER®

TRABECULAR LASER MELTED TITANIUM




permedica
ORTHOPAEDICS

TRASER®

Trabecular Titanium
by Selective Laser Melting



TRASER[®] Trabecular lattice




“Mimicking the nature of bone to provide bone ingrowth and osseointegration for long-term implant stability”

“Highly porous trabecular network characterized by open fully interconnected irregular pores to promote fast bone ingrowth”




Porosity

Void volume - Porosity: 70%.
Random irregular shaped pores.
Complete permeability.
Pore interconnectivity: 100%.
Pore size range: 100-2000 μm .



Isotropy


Random pore sizes and random irregular trabeculae orientation and distribution along the 3-dimensions.
Isotropic uniform highly-porous lattice.



“Highly osteoconductive properties”

Superficial topography

A glimpse on the topography of the trabeculae reveals a uniform microroughened surface (Ra, 2 μm) which optimize osteointegration (bone on-growth).



Surface Roughness

It has been demonstrated that the bone response is influenced by the implant surface topography at the micrometer level. Surface roughness affects the physiological processes of bone tissue synthesis and osseointegration. Osteoblasts response in terms of cell adhesion, proliferation and osteosynthesis local factor production, as well as implant-to-bone direct contact have been proven to be increased by moderately-roughened surfaces (Ra, 1 to 6 μm) [1-11].

Fully interconnected open irregular-shaped pores

Cell/tissue ingrowth behavior depends on the pore size, porosity and pore interconnectivity in different ways. Pores have to be open and fully interconnected to allow cell migration and to permit vascularization for cell viability and proliferation. Formation of new bone into an interconnected network of pores provide mechanical interlocking between implant and the surrounding bone [9,10].

High porosity


Higher porosity means higher bone ingrowth. A large consensus suggests a minimum porosity of 60% to maximize bone ingrowth while maintaining at the same time adequate mechanical properties of the porous structure [12,13,14].

Mean pore size of 500 micron

New bone formation is promoted by a pore size range of 150-500 μm . Smaller pore sizes leads to fibrous tissue ingrowth and prevent the mineralization of the osteoid matrix. Greater pore sizes lead to a slower bone ingrowth and incomplete pore filling [15-20].


Mechanical Properties

Compressive strength of TRASER® lattice: 38 Mpa




High friction coefficient of TRASER® lattice with cancellous bone to improve implant primary fixation [21-24].


Elastic modulus [Gpa]: 1.4 [1.1-1.7]



Coefficient of friction of TRASER® lattice with cancellous bone: 2.1



TRASER® is not a coating




Traser® is not a porous coating. The solid and the porous portions of the cup are built up in a one continuous process, thus creating a single piece without interface layers between different portions and, consequently, without shear forces.



Jump System *TRASER*® Cup



SLM



Selective Laser Melting

Additive Manufacturing

Additive manufacturing (or 3D printing) is a technology that allows to produce three-dimensional solid objects from a digital model using as a raw material metal powder.

Does not replace completely but rather complements the traditional machining processes.

The objects are created by adding the material in successive layers one over the other up to completion of the object.

This differs from traditional machining process with machine tools, in which the object is obtained by subtracting material from the raw material.

Selective Laser Melting is an additive manufacturing technology that selectively melts and sinterizes by means of the thermal energy from a laser beam specific portions of titanium powder layer to create 3D solid parts

Selective Laser Melting process take place in an inert atmosphere (Argon) in order to avoid any titanium powder oxidation.

TRASER[®] Raw material Titanium powder - Ti6Al4V (ISO 5832/3)


CHEMICAL ANALYSIS

Element	Result (Wt %)
Aluminium	6.13
Carbon	0.01
Iron	0.21
Hydrogen	0.002
Nitrogen	0.034
Oxigen	0.09
Titanium	Balance
Vanadium	3.82


Production flow




3D - CAD Model



Titanium Powder



power diameter range: 16 - 76 µm



SLM



Green State



Thermal Treatment Sandblasting



Traditional milling and turning processes

Customization

A clinical case with custom made components

HISTORY

During the first consultancy in November 2015 in the hospital, 32y old male patient; dwarfism; 1.08m/30 kg of bodyweight.

In the anamnesis continuous daily pain in the right hip with restricted range of movement. Important limitation in QDL (quality of daily living). Pain killers dependent.

CLINICAL EXAM

Walk : limping right side; Duchenne/Trendelenburg negative; ROM : F/E 30-10-0, ER/IR 20-0-10, ABD/ADD 10-0-20.

Radiologically : secondary osteoarthritis of the right hip (Crowe II, Paprovsky IIb).

SURGERY

April 17th, 2016, total hip replacement is performed, using an allograft out of the femoral head, fixed with several screws to build up the lateral acetabular rim. Then a TRASER® Trabecular Titanium custom made cup (Permedica Orthopaedics, Merate, Italy) is inserted.

On the femoral side, a customized selective laser melted femoral stem is implanted; reposition with a metal femoral head.


The final result is stable with a perfect covering of the cup implant and central seating of the TRASER® custom made cup.

FOLLOW UP

4 months after surgery perfect walk. No crutches since the second month postoperatively; no limping.

Joint function : F/E 120-0-20, ER/IR in 90° of flexion 40-0-30. Radiologically the implant is stable; bone graft is stable, the cup is integrated.

“Freedom of design to meet patient/surgeon requirements.”



Actual Size!



Pre-Op



Post-Op



6-Months

References

1. Feighan JE, Goldberg VM, Davy D, Stevenson S. The influence of surface-blasting on the incorporation of titanium-alloy implants in a rabbit intramedullary model. *JBS Am.* 1995;77:1380-95.
2. Wennerberg A, Albrektsson T. Effects of titanium surface topography on bone integration: a systematic review. *Clin Oral Implants Res.* 2009 Sep;20 Suppl 4:172-84.
3. Lincks J, Boyan BD, Blanchard CR, Lohmann CH, Liu Y, Cochran DL, et al. Response of MG63 osteoblast-like cells to titanium and titanium alloy is dependent on surface roughness and composition. *Biomaterials* 1998;19:2219-32.
4. Herrero-Climent M, Lazaro P, Vicente Rios J, Lluch S, Marques M, Guillem-Marti J, et al. Influence of acid-etching after grit-blasted on osseointegration of titanium dental implants: in vitro and in vivo studies. *J Mater Sci Mater Med* 2013;24:2047-55.
5. Ivanoff CJ, Hallgren C, Widmark G, Sennerby L, Wennerberg A. Histologic evaluation of the bone integration of TiO₂ blasted and turned titanium microimplants in humans. *Clin Oral Implants Res* 2001;12:128-34.
6. Gotfredsen K, Wennerberg A, Johansson C, Skovgaard LT, Hjorting-Hansen E. Anchorage of TiO₂-blasted, HA-coated, and machined implants: an experimental study with rabbits. *J Biomed Mater Res* 1995;29:1223-31.
7. Wennerberg A, Albrektsson T, Johansson C, Andersson B. Experimental study of turned and grit-blasted screw-shaped implants with special emphasis on effects of blasting material and surface topography. *Biomaterials* 1996;17:15-22.
8. De Wild M, Schumacher R, Mayer K, Schkommodau E, Thoma D, Bredell M, Kruse Gujer A, Grätz KW, Weber FE. Bone regeneration by the osteoconductivity of porous titanium implants manufactured by selective laser melting: a histological and micro computed tomography study in the rabbit. *Tissue Eng Part A.* 2013 Dec;19(23-24):2645-54.
9. Mour M, Das D, Winkler T, Hoenig E, Mielke G, Morlock M, Schilling A. Advances in porous biomaterials for dental and orthopaedic applications. *Materials.* 2010;3:2947-74.
10. Nouri A, Hodgson PD, Wen C. Biomimetic porous titanium scaffolds for orthopedic and dental applications. In Mukerjee A, Ed. *Biomimetics Learning from Nature.* InTech 2010;21:415–50.
11. Hacking SA, Bobyn JD, Tanzer M, Krygier JJ. The osseous response to corundum blasted implant surfaces in a canine hip model. *Clin Orthop Relat Res.* 1999 Jul;(364):240-53.
12. Bragdon CR, Jasty M, Greene M, Rubash HE, Harris WH. Biologic fixation of total hip implants. Insights gained from a series of canine studies. *JBS Am.* 2004;86-A Suppl 2:105-117.
13. Karageorgiou V, Kaplan D. Porosity of 3D biomaterial scaffolds and osteogenesis. *Biomaterials.* 2005 Sep;26(27):5474-91.
14. Cheng A, Humayun A, Cohen DJ, Boyan BD, Schwartz Z. Additively manufactured 3D porous Ti-6Al-4V constructs mimic trabecular bone structure and regulate osteoblast proliferation, differentiation and local factor production in a porosity and surface roughness dependent manner. *Biofabrication.* 2014 Oct 7;6(4):045007.
15. Warnke PH, Douglas T, Wollny P, Sherry E, Steiner M, Galonska S, Becker ST, Springer IN, Wiltfang J, Sivananthan S. Rapid prototyping: porous titanium alloy scaffolds produced by selective laser melting for bone tissue engineering. *Tissue Eng Part C Methods.* 2009 Jun;15(2):115-24.
16. Simske SJ, Ayers RA, Bateman TA. Porous materials for bone engineering. *Mater Sci Forum.* 1997;250:151–182.
17. Bobyn JD, Pilliar RM, Cameron HU, Weatherly GC. The optimum pore size for the fixation of porous-surfaced metal implants by the ingrowth of bone. *Clin Orthop Relat Res.* 1980 Jul-Aug;(150):263-70.
18. Miao X, Sun D. Graded/Gradient Porous Biomaterials. *Materials.* 2010;3:1-22.
19. Bansiddhi A, Sargeant TD, Stupp SI, Dunand DC. Porous NiTi for bone implants: A review. *Acta Biomater.* 2008 Jul;4(4):773–782.
20. Frosch KH, Barvencik F, Viereck V, Lohmann CH, Dresing K, Breme J, Brunner E, Stürmer KM. Growth behavior, matrix production, and gene expression of human osteoblasts in defined cylindrical titanium channels. *J Biomed Mater Res A.* 2004 Feb 1;68(2):325-34.
21. Zhang Y, Ahn PB, Fitzpatrick DC, Heiner AD, Poggie RA, Brown TD. Interfacial frictional behavior: Cancellous bone, cortical bone, and a novel porous tantalum biomaterial. *Journal of Musculoskeletal Research.* 1999;3(4):245-251.
22. Shirazi-Adl A, Dammak M, Paiement G. Experimental determination of friction characteristics at the trabecular bone/porous-coated metal interface in cementless implants, *J Biomed Mater Res.* 1993;27:167-175.
23. Rancourt D, Shirazi-Adl A, Drouin G, Paiement G (1990) Friction properties of the interface between porous-surfaced metals and tibial cancellous bone. *J Biomed Mater Res* 24:1503–1519
24. Janssen D, Zwartelé RE, Doets HC, Verdonchot N. Computational assessment of press-fit acetabular implant fixation: the effect of implant design, interference fit, bone quality, and frictional properties. *Proc Inst Mech Eng H.* 2010;224(1):67-75.
25. Internal technical file Permedica.

CHALLENGING EXCELLENCE IN TECHNOLOGY



permedica
ORTHOPAEDICS

permedica s.p.a.- Cap. Soc. 1.300.000 €
via Como, 38 - 23807 Merate (LC) ITALY
tel.: +39 039 95.14.811 - fax +39 039 99.03.078
www.permedica.com - info@permedica.com