Medtronic

Pacific™ Plus

PTA Catheter
PTA катетър
Katéтр pro PTA
PTA-kateter
PTA-Katheter
Καθετήρας PTA (διαδερμικής διαυλικής αγγειοπλαστικής)
Catéter para ATP
PTA kateeter
PTA-katetri
Cathéter pour angioplastie transluminale percutanée (ATP)
Kateter za PTA
PTA katéter
Kateter PTA
Catetere per PTA
ТТА катетері
PTA 카테터
PTA kateteris
PTA katetrs
PTA катетер (катетер за перкутана транслуминална ангиопластика)
PTA-katheter
PTA-kateter
Cewnik do PTA
Cateter para PTA
Cateter ATP
Катетер баллонный для ЧТА
PTA katéter
Kateter PTA
PTA kateter
PTA-kateter
PTA Kateteri
Катетер для ЧТА
Ống thông PTA

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Petunjuk Penggunaan • Istruzioni per l'uso • Қолдану жөніндегі нұсқаулар • 사용 지침 • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з використання • Hướng dẫn sử dụng

CE
2797
UK
CA
0086

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Contents

1 Description	12
2 Intended purpose	12
2.1 Intended patient population	12
2.2 Indications for use	12
2.3 Contraindications	12
2.4 Clinical benefits	12
2.5 Intended users	13
2.6 Device performance characteristics	13
3 Warnings	13
4 Precautions	13
5 Adverse effects	14
6 Instructions for use	14
6.1 Inflation device preparation	14
6.2 Catheter selection	14
6.3 Preparation technique	15
6.4 To displace air	15
6.5 Insertion technique	15
6.6 Balloon inflation	15
6.7 Balloon deflation and removal	16
7 Device disposal	16
8 How supplied	16
9 Storage	16
10 Disclaimer of warranty	16

1 Description

The Pacific Plus PTA catheter is an Over the Wire (OTW) peripheral balloon catheter, specifically designed for percutaneous transluminal angioplasty (PTA) in stenosed vessel segments. The OTW catheter is used to guide the balloon to the stenosed vessel segment. The balloon is then inflated to dilate the vessel.

The catheter is a coaxial dual lumen device. The lumen marked “WIRE” is the central lumen of the catheter, which terminates at the distal tip. This lumen is used to pass the catheter over a guidewire with a maximum outer diameter of 0.018 in (0.46 mm). The lumen marked “BALLOON” is the balloon inflation lumen, which is used to inflate and deflate the dilatation balloon with a mixture of contrast medium and saline solution.

Consult the Compliance Chart included in the catheter packaging for the diameters of the balloons at given pressures. The Pacific Plus PTA catheter is available in different balloon sizes and is provided with a hydrophilic coating. The balloon length is specified on the package label.

The Pacific Plus PTA catheter does not contain natural rubber latex. However, during the manufacturing, assembly, or packaging process, it may have incidental contact with materials, tools, or equipment containing or coming in contact with latex.

2 Intended purpose

The intended purpose of the Pacific Plus PTA catheter is to restore lumen patency and blood flow to peripheral arteries, or through native arteriovenous dialysis fistulae.

The Pacific Plus PTA catheter is not for use in direct contact with the heart or the central circulatory system.

2.1 Intended patient population

The Pacific Plus PTA catheter is intended to improve and maintain luminal diameter in patients at risk for closure or restenosis following percutaneous transluminal angioplasty in peripheral arteries or who have obstructive lesions in native arteriovenous dialysis fistulae.

This device has not been tested in a pediatric population.

2.2 Indications for use

The Pacific Plus PTA catheter is intended to dilate stenoses in the iliac, femoral, iliofemoral, popliteal, and infrapopliteal arteries, and for the treatment of obstructive lesions of native arteriovenous dialysis fistulae.

2.3 Contraindications

The catheter is contraindicated for use in coronary arteries or the neurovasculature. It is also contraindicated when the guidewire is unable to cross the target lesion.

2.4 Clinical benefits

The clinical benefits of treating stenoses in peripheral arteries or obstructive lesions in native arteriovenous dialysis fistulae with the Pacific Plus PTA catheter include the following:

- Restored lumen patency and blood flow to resolve arterial disease-related symptoms
- Improved quality of life
- Increased mobility

The Pacific Plus PTA catheter has been found to be safe and effective with clinical benefits outweighing the risks in the intended population.

Clinical data demonstrating these benefits are shown in *Table 1* and *Table 2* below.¹

Table 1. Functional endpoints at 1 year

	Rate (n/N) [95% CI]
Functional status ^a improved from baseline	22.2% (4/18) [6.4%, 47.6%]
Functional status preserved from baseline	38.9% (7/18) [17.3%, 64.3%]
Ambulation status ^b improved from baseline	11.8% (2/17) [1.5%, 36.4%]
Ambulation status preserved from baseline	70.6% (12/17) [44.0%, 89.7%]

¹ Data are sourced from the Society of Vascular Surgery Vascular Quality Initiative (VQI) registry of vascular procedures performed between July 2018 and August 2021 by 540 physicians at 188 participating centers in the United States of America. A total of 70 patients and 83 arteries were treated by percutaneous transluminal angioplasty (PTA) with the Pacific Plus PTA catheter being the only device used for the treatment of peripheral artery disease. One-year follow-up data were available for 22 patients, and there were data pertinent to the clinical benefits available in 17 or 18 of these patients, depending on the endpoint shown above.

Table 1. Functional endpoints at 1 year (continued)

	Rate (n/N) [95% CI]
Rutherford classification improved ≥ 1 classification from baseline	77.8% (14/18) [52.4%, 93.6%]
Rutherford classification preserved from baseline	11.1% (2/18) [1.4%, 34.7%]

^a Functional status categories (from worst to best) are assisted care, self-care, light work, and full function.

^b Ambulation status categories (from worst to best) are bed bound, wheelchair only, ambulation with assistance, ambulation with prosthesis, and full ambulation.

Table 2. Kaplan-Meier estimate^a of primary patency at 1 year

	Estimate $\pm$ SE (N at risk) [95% CI]
Primary patency of treated artery	89.8% $\pm$ 3.9% (N=83) [82.2%, 97.4%]

^a Kaplan-Meier estimate of primary patency at one year from an initial number-at-risk of 83 arteries. Primary patency is defined as freedom from flow-limiting stenosis or occlusion without surgical or interventional revascularization in the treated lesion following the initial PTA procedure.

2.5 Intended users

Pacific Plus PTA catheter percutaneous transluminal angioplasty procedures should only be performed by physicians who have experience with interventional techniques in the vascular system.

The device should only be used within a clinical setting in a sterile environment.

2.6 Device performance characteristics

The Pacific Plus PTA catheter treats occlusions or lesions in peripheral arteries or in native dialysis fistulae. The balloon is inflated to its predetermined diameter and exerts a constant, gentle outward force to establish patency

3 Warnings

- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation.
- Do not manipulate the PTA balloon in an inflated state. The position of the PTA balloon may only be changed with the guidewire in place.
- If resistance occurs during manipulation, determine the cause using fluoroscopy, roadmapping, or DSA before moving the catheter backward or forward.
- The guidewire may under no circumstances be moved during inflation of the balloon.
- The balloon must be completely deflated before removing the catheter from the vascular system.
- Do not exceed the rated burst pressure (RBP). The rated burst pressure is based on the results of in vitro testing. Use of a pressure monitoring device is recommended to prevent over pressurization. Inflation in excess of the rated burst pressure may cause the balloon to rupture.
- Use a 1:1 or 1:3 mixture of contrast medium and saline solution to fill the balloon. (It has been shown that a 1:3 ratio yields faster balloon inflation/ deflation times.) Never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to the Use By date noted on the package.
- Do not expose to organic solvents (such as alcohol), ionizing radiation, or ultraviolet light.
- The length of the introducer sheath should be appropriate to the support needed, and should be compatible with the PTA catheter characteristics listed on the package label.

4 Precautions

- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- Use an introducer sheath with a French compatibility equal to or larger than the one reported in the label. Using an inappropriate introducer sheath may affect the advancement or retrieval of the catheter.
- Appropriate drug therapy (anticoagulant, vasodilator, etc.) should be administered to the patient according to standard protocols for PTA before insertion of the balloon dilatation catheter.

Note: Use appropriate anticoagulation to prevent or reduce device thrombosis, as this device has not been tested in the absence of anticoagulation.

- Exercise care during handling to avoid possible damage to the catheter. Do not use a catheter that has been damaged.

- To minimize the possible introduction of air into the system, before proceeding pay careful attention to the maintenance of tight catheter connections throughout aspiration and flushing of the system.
- Exercise caution when using the Pacific Plus PTA catheter for procedures with calcified lesions due to the abrasive nature of these lesions.
- Allergic reactions to contrast medium should be identified before treatment.
- The general technical requirements for catheter insertion must be observed at all times. This includes flushing the components with sterile, isotonic saline solution prior to use and the usual prophylactic, systemic heparinization.
- Catheter applications vary. Select the technique to be used based on the patient's condition and the experience of the interventionalist.
- The maximum diameter of the guidewire must not exceed 0.018 in (0.46 mm). Certain coatings on guidewires may lead to increased friction between the guidewire and the catheter.

5 Adverse effects

Complications associated with the use of the Pacific Plus PTA catheter are similar to those associated with standard PTA procedures. Possible complications may include, but are not limited to:

Puncture related:

- Arteriovenous fistula
- Local hematoma
- Local hemorrhage
- Local infections
- Local or distal thromboembolic episodes
- Pseudoaneurysm
- Thrombosis

Dilatation related:

- Acute reocclusion necessitating surgical intervention
- Dissection in the dilated artery wall or fistula
- Perforation of the artery wall
- Prolonged spasms
- Restenosis of the dilated artery
- Total occlusion of the peripheral artery

Angiography related:

- Allergic reaction to contrast medium
- Arrhythmias
- Death
- Drug reactions
- Endocarditis
- Hypotension
- Pain and tenderness
- Sepsis/infection
- Short-term hemodynamic deterioration
- Systemic embolization

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

6 Instructions for use

6.1 Inflation device preparation

Prepare the inflation device according to the manufacturer's instructions.

6.2 Catheter selection

1. The nominal balloon size must be chosen equal to or smaller than the inner diameter of the artery proximal and distal to the lesion.

2. If the stenosis cannot be crossed with the desired catheter, use a smaller catheter to predilate the lesion and facilitate the passage of a more appropriately-sized catheter.

6.3 Preparation technique

Caution: Avoid manipulating the catheter while removing it from the packaging and flushing the guidewire lumen.

1. Prior to use, carefully examine the unit to verify that the catheter or sterile package has not been damaged in shipment.
2. The catheter is packaged in a protective coil. Gently remove it from the hoop.
3. After moistening the balloon protection sheath with physiological saline solution, carefully remove the balloon protection sheath without damaging the balloon.

Caution: Do not use with Lipiodol or Ethiodol contrast media, or other such contrast media that incorporate the components of these agents.

6.4 To displace air

1. Flush out the guidewire lumen by connecting a syringe filled with saline solution. Start flushing until liquid rinses out of the guidewire lumen at the distal tip of the catheter.
2. Connect a syringe containing contrast medium to the luer connector located at the proximal hub of the catheter and apply negative pressure for approximately 15 seconds until no bubbles appear in the contrast medium solution. Return to atmospheric pressure, allowing contrast medium flow into the catheter lumen. Remove the syringe, leaving a meniscus of contrast medium in the hub of the balloon lumen.
3. Remove air bubbles from the inflation device following the manufacturer's instructions.
4. Using the stopcock, attach the inflation device to the catheter. Do not allow air to enter the system.

Note: The Pacific Plus PTA catheter is provided with a hydrophilic coating. Prior to inserting the catheter, activate the coating by immersing the catheter in normal saline for approximately 5 seconds or wiping down the catheter shaft with a saline-saturated gauze.

Caution: Do not wipe down the catheter surface with a dry gauze.

6.5 Insertion technique

Note: Access to the vessel must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention.

1. If the guidewire has not already been positioned across the lesion, under fluoroscopy insert a 0.018 in (0.46 mm) or 0.014 in (0.36 mm) guidewire across the lesion following standard PTA techniques.
2. Insert the proximal end of the guidewire into the distal tip of the catheter.
3. Carefully insert the catheter through the introducer sheath and, under fluoroscopic guidance, carefully advance it through the peripheral vasculature following the guidewire toward the lesion.

Caution: Do not advance the guidewire or the catheter if resistance is met without first determining the cause of resistance and taking remedial action.

4. Carefully position the balloon within the lesion using the 2 radiopaque markers located under the balloon as reference points for precise placement across the target lesion. For balloons ≥ 200 mm, there are 2 additional radiopaque markers at the center of the balloon.
5. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to predilate the lesion and facilitate passage of a more appropriately-sized dilatation catheter.

Note: To avoid kinking, advance the catheter slowly in small increments until the proximal end of the guidewire emerges from the catheter.

6.6 Balloon inflation

1. Make sure all air bubbles are removed from the catheter and the inflation device before inflating the balloon.
2. Open the stopcock on the inflation device. Inflate the balloon gradually to dilate the lesion to the calculated diameter in accordance with the Compliance Chart.
3. To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.

Caution: Do not inflate the balloon more than 10 times at the RBP.

4. Apply a constant pressure for approximately 30 seconds.

Caution: Do not exceed the rated burst pressure (RBP).

Note: The subject device was tested on the bench with the Protégé™ EverFlex™ (Medtronic) self-expanding peripheral stent system. All stents should be deployed in accordance with the manufacturer's indications and instructions for use.

6.7 Balloon deflation and removal

Caution: Larger models of the Pacific Plus PTA catheter may exhibit slower deflation times, particularly on long catheter shafts.

1. Deflate the balloon in accordance with standard PTA procedures. Apply negative pressure to the balloon for at least 60-120 seconds. Ensure the balloon is fully deflated under fluoroscopy before carefully removing the catheter from the vessel.
2. Carefully remove the deflated balloon from the introducer while maintaining negative pressure.
3. After use, this product may be a potential biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable local, regional, and national laws and regulations.

7 Device disposal

Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

8 How supplied

The Pacific Plus PTA catheter is supplied sterile and intended for single use only. It is sterilized by ethylene oxide gas and it will remain sterile as long as the packaging remains unopened and undamaged. Use product prior to Use By date. Do not use if the packaging is damaged.

Caution: Do not use if the inner package is open or damaged.

9 Storage

Store the device in a dry area away from sunlight.

Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used before the Use-by date on the package label.

10 Disclaimer of warranty

Note: This disclaimer of warranty does not apply in Australia.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
**Medtronic, Inc.**
Medtronic B.V.
Medtronic Limited