# RealPCR\* ASFV DNA Test

## English Version

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**REF** 99-56020

 **Version**  
06-56020-04

*This insert is not for a product available to U.S. or Canadian customers.  
Please refer to the RealPCR insert webpage and select "United States  
or Canada" to view the correct product insert.*

Test With Confidence™

**IDEXX**

MENU

PRINT

English version

# RealPCR\* ASFV DNA Test

For veterinary use only

## Name and Intended Use

The RealPCR ASFV DNA Test is used for the detection of African swine fever virus (ASFV) DNA extracted from serum, plasma, EDTA-blood, oral fluids, tissue (lung, spleen, kidney, bone marrow, lymph nodes and tonsils) and swab samples from swine, including wild boars. Samples can be tested in pools of up to 20. Pools containing a single weak sample (for example, Ct > 32) may yield a negative result due to the dilution effect of pooling.

## General Information

African swine fever virus is the causative agent of haemorrhagic African swine fever. ASFV is a highly contagious agent, that with subclinical symptoms can persistently infect warthogs, bushpigs as well as *Onnithodoros* soft ticks. As a reportable disease to the World Health Organization for Animal Health (OIE), ASFV has a high morbidity and mortality in domestic swine. Controlling ASFV outbreaks relies heavily on early detection of the virus. The RealPCR ASFV DNA Test provides a highly sensitive and specific means to detect viral DNA.

The IDEXX RealPCR system is a modular format in which disease-specific target mixes are paired with a standardized DNA or RNA master mix and a single pooled positive control. Reagents are individually packaged and sold separately to allow for flexible reagent handling.

The RealPCR ASFV DNA Mix (ASFV DNA Mix) contains primers and probes for the detection of ASFV DNA when amplified with RealPCR DNA Master Mix (DNA MMx). The internal control for the test is based on the detection of a genomic DNA sequence in the sample that is conserved in swine. This genomic target is referred to as the internal sample control (ISC) in this protocol. Detection of endogenous DNA in the test sample controls for sample addition, extraction and amplification. Primers and probe for detection of the internal sample control are included in the ASFV DNA Mix. An optional internal positive control, the RealPCR Internal Positive Control (IPC  $\geq$  v1.2), is also available and should be used when endogenous host DNA is at low levels or unlikely to be present after extraction (such as feed samples). The IPC contains a synthetic version of the swine ISC DNA target and is therefore compatible with the ASFV DNA Mix. Refer to the RealPCR Internal Positive Control ([REF](#) 99-56330) product insert for guidance.

## Materials and Storage

Identification/ General Information	Cap color	Quantity	Storage		Freeze/Thaw cycles
		100 tests	At receipt	After reconstitution	
<b>RealPCR ASFV DNA Mix (ASFV DNA Mix), dried</b>	Green	1 x 1.0 mL	-25 to 8°C	-25 to -15°C	≤6
<p><b>REF</b> 99-56020</p> <p>Reconstitute to 1 mL in PCR Grade Water. Store the ASFV DNA Mix in the dark. The expiration date on the vial is valid for either the dry or reconstituted form. The ASFV DNA Mix tube label indicates the PC version that is compatible with that target mix. For example, PC ≥ v1.4 means the target mix can be used with PC versions greater or equal to 1.4.</p>					
<b>RealPCR DNA Master Mix (DNA MMx)</b>	Purple	1 x 1.0 mL	-25 to -15°C (Long-term)	N/A	≤6
<p><b>REF</b> 99-56250</p> <p>Concentrated master mix that includes hot-start polymerase for use with DNA target mixes in the IDEXX RealPCR system. The DNA MMx is more viscous than most master mixes– see the Test Procedure section for handling recommendations. A reference dye (ROX) has been added for normalizing volume inaccuracies. Protect the DNA MMx from light.</p>					
<b>RealPCR Positive Control, dried (PC)</b>	Blue	1 x 500 µL	-25 to 8°C	-25 to -15°C	≤6
<p><b>REF</b> 99-56310</p> <p>Reconstitute to 500 µL in PCR Grade Water. The PC contains all IDEXX RealPCR and ISC targets (including the target for ASFV) and is intended for use with all IDEXX RealPCR target mixes. The expiration date on the vial is valid for either the dry or reconstituted form. The PC has a version number (such as v1.3). When new target mixes are developed for the RealPCR product line, the target sequences are added to the PC and the PC version number increases (for example, v1.3 would be updated to v1.4).</p> <p>The PC includes the IDEXX Signature (a unique oligonucleotide sequence). Presence of the IDEXX Signature, in the work environment, indicates PC contamination. For laboratories that want to monitor for PC contamination, the IDEXX Signature can be detected using the RealPCR PC Tracker DNA Mix and RealPCR DNA MMx.</p>					
<b>RealPCR PCR Grade Water</b>	Clear	2 x 1.0 mL	-25 to 8°C		N/A
<p><b>REF</b> 99-56350</p> <p>PCR Grade Water has been qualified for reverse transcription-PCR (RT-PCR) use. It is used for the reconstitution of RealPCR reagents. It is also used as the PCR negative control for each test run. Do not transfer PCR Grade Water vials between PCR work areas. Separate vials of water are needed for each area to avoid contamination risk.</p>					

**Note:** See table at the end of the insert for a description of symbols used on the insert and labels.

## Materials Required but Not Provided

- Commercial DNA extraction kit
- Optional—Centrifuge with a rotor and adapters for multi-well plates
- Micro-centrifuge for 2 mL microtubes capable of reaching 1500 – 3000 x g
- Appropriate personal protective equipment (e.g., gloves, lab coat)
- Nuclease-free, aerosol-resistant pipette tips
- Sterile microtubes for preparation of PCR Mix
- Pipettes (5–1000  $\mu$ L); dedicated pipettes for preparation of PCR Mix
- 96- or 384-well format PCR plates and optical adhesive film/plate covers

- Real-time PCR instrument

(Applied Biosystems® 7500, Applied Biosystems® 7500 Fast System (Standard and Fast Mode), Applied Biosystems® ViiA™ 7, Applied Biosystems QuantStudio 5, Agilent Mx3000P™, Agilent Mx3005P™, Agilent AriaMx, Bio-Rad CFX96 Touch™, Bio Molecular Systems Mic qPCR Cycler, QIAGEN Rotor-Gene (72-Well Rotor only), Roche LightCycler® 480 or equivalent).

**Note-** the Roche LC480 instrument requires additional calibration and software settings.

IDEXX Technical Service can provide guidance for use of the instruments mentioned above with RealPCR reagents.

## Laboratory Practices and Warnings

- Do not use reagents past expiration date.
- The entire procedure must be performed under nuclease-free conditions.
- Wear powder-free gloves when working with the reagents and nucleic acids.
- To avoid cross-contamination, use nuclease-free, aerosol-resistant pipette tips for all pipetting, and physically separate the workplaces for nucleic acid extraction/handling, PCR setup and PCR.

## Reconstitution of Dried Components

Reconstitute the ASFV DNA Mix and Positive Control by pipetting PCR Grade Water to the volume indicated on the component label. Allow to sit at 18 to 26°C for at least 10 minutes; mix and microcentrifuge briefly prior to use. Once the ASFV DNA Mix and the Positive Control are reconstituted, aliquot as appropriate and store the solutions frozen. When handling frozen components, thaw at 18 to 26°C for approximately 15 to 30 minutes, mix gently and then microcentrifuge briefly (~1,500 – 3,000  $\times$  g).

## DNA Extraction

The ASFV DNA Mix has been validated using the commercial extraction methods listed below. Other extraction or lysis methods may also be used once validated by the laboratory.

RealPCR DNA/RNA Magnetic Bead Kit (IDEXX)  
RealPCR DNA/RNA Spin Column Kit (IDEXX)  
MagMAX Pure- 96 Viral RNA Isolation Kit (Thermo)  
QIAamp DNA Mini Kit (Qiagen)  
High Pure PCR Template Prep Kit (Roche)

Store the purified DNA at  $\leq -15^{\circ}\text{C}$  if testing is not performed immediately after DNA extraction. It is recommended that a negative extraction control (mock sample) is included as a sample.

## Test Procedure

- 1 Preparation of the PCR Mix.
  - Mix the thawed DNA MMx by inversion or gentle vortex.
  - The DNA MMx is a viscous solution; always pipette it slowly.
  - To prepare the PCR Mix add 10  $\mu$ L ASFV DNA Mix and 10  $\mu$ L DNA MMx for each reaction.
  - When preparing the PCR Mix, first pipette ASFV DNA Mix into the tube and then add the DNA MMx. Pipette up and down a few times to rinse the MMx pipette tip.
  - Gently vortex the solution to ensure the components are mixed well.
  - Pipette the PCR Mix slowly into the PCR plate.

The PCR Mix can be stored at 2 to 8°C for up to 24 hrs or –25 to –15°C for up to 2 weeks. Protect from light.

- 2 Pipette 20  $\mu$ L of the PCR Mix into the required wells of the multiwell plate.
- 3 Add 5  $\mu$ L of sample DNA to each well. The final reaction volume is 25  $\mu$ L.
- 4 Include the Positive Control (5  $\mu$ L) and PCR negative control (5  $\mu$ L PCR Grade Water) for each test run.
- 5 Cover the plate and briefly spin the plate, if necessary, to settle contents and remove air bubbles.
- 6 Set up thermal cycler with IDEXX RealPCR Standard DNA/RNA Cycling Program.

### Settings for Reporter and Quencher

Target	Reporter	Quencher
ASFV	FAM™	BHQ® (none)
Internal Control (ISC/IPC)	HEX™ (VIC)	BHQ (none)
Passive Reference	ROX™	N/A

### RealPCR\* Standard DNA/RNA Cycling Program

	Temperature	Time	Cycles
Reverse transcription (RT)	50°C	15 min.	1
Denaturation	95°C	1 min.	1
Amplification**	95°C	15 sec.	45
	60°C	30 sec.	

\*\*Ensure the instrument is set to record fluorescence following the 60°C amplification step.

**Note-** for DNA targets, the cycling protocol can be run without the RNA "RT" step. However it is recommended that the "RT" step is used routinely so that RNA tests can easily be incorporated into the work-flow.

---

## 7 Analyze data

When setting up the instrument software, assign a unique identifier to each target and internal control. For example, when target A and B are on the same plate, the A wells must be analyzed independently from the B wells. Refer to specific instrument's user manual for guidance on how to analyze data. To set the threshold, use the Auto Ct setting.

- Agilent Mx3000P and Mx3005P– ensure the background-based threshold fluorescence method is being used for analysis.
- QIAGEN Rotor-Gene instrument– manually adjust the threshold line above background in the linear phase of exponential amplification. This is best done on log view plots and should be repeated for each reporter in the target mix.
- Applied Biosystems instruments– In certain situations the auto threshold setting does not produce satisfactory results. In these situations, a manual adjustment of the threshold is necessary to determine Ct values. This is best done on log view plots and should be repeated for each reporter in the target mix.

### Validity Criteria

	<u>FAM Ct value</u>	<u>HEX (VIC) Ct value</u>
Positive Control	<38	<38
PCR Negative Control	No Signal	No Signal

---

## 8 Interpretation of results

<u>Sample Result</u>	<u>FAM Signal</u>	<u>HEX (VIC) Signal</u>	<u>Other characteristics</u>
ASFV DNA detected	Yes	Yes/No	A positive Ct value and characteristic amplification curve in comparison to the PCR negative control. An internal control amplification curve in the HEX (VIC) channel is expected; some strong ASFV-positive DNA samples may result in a negative internal control result. <sup>†</sup>
ASFV DNA not detected	No	Yes	Amplification curve in the HEX (VIC) internal control channel.
Invalid <sup>‡</sup>	No	No	Absence of an amplification curve in the FAM and HEX (VIC) channels indicates an invalid result.

<sup>†</sup>The target mix is optimized for the detection of ASFV DNA; a strong positive DNA sample may out compete the detection of the internal control.

<sup>‡</sup>An invalid sample can be an indication of failed sample addition, extraction and/or PCR. It is recommended that the DNA be diluted five-fold into PCR grade water and retested; include the undiluted DNA as a sample. If the test is still not valid a new extraction is recommended.

**For technical assistance:**

IDEXX USA Tel: +1 800 548 9997 or +1 207 556 4895

IDEXX Europe Tel: +800 727 43399

Contact your IDEXX area manager or distributor or visit our website: [idexx.com/contact/pd](https://idexx.com/contact/pd)

\*IDEXX, RealPCR and Test With Confidence are trademarks or registered trademarks of IDEXX Laboratories or its affiliates in the United States and/or other countries. All other product and company names and logos are trademarks of their respective holders.

Dye compounds in this product are sold under license from Biosearch Technologies, Inc. and protected by U.S. and world-wide patents either issued or in application. The license covers veterinary applications. Not for human IVD use.

Patent information: [idexx.com/patents](https://idexx.com/patents)

©2022 IDEXX Laboratories, Inc. All rights reserved.

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IDEXX Montpellier SAS
IDEXX Europe B.V.