IRISpec CA/CB/CC

REF 800-7702

TABLE OF CONTENTS

English	2
Français (France) (fr-fr)	6
Deutsche (de-de)	10
Italiano (it-it)	14
Español (es-es)	18
Português (Portugal) (pt-pt)	22
Dansk (da-dk)	26
Svenska (sv-se)	30
Norsk (nb-no)	34
Ελληνικά (el-gr)	38
日本語 (ja-jp)	42
中文 (中国) (zh-cn)	46
Magyar (hu-hu)	50
Polski (pl-pl)	54
Čeština (cs-cz)	58
Slovák (sk-sk)	62
한국어 (ko-kr)	66
Türkçe (tr-tr)	70
Русский (ru-ru)	74
Hrvatski (hr-hr)	79
Български (bg-bg)	83
中文 (台灣) (zh-tw)	87
Română (ro-ro)	91
Srpski (sr-rs)	95
Українська (uk-ua)	99
Português (Brasil) (pt-br)	103
Tiếng Việt (vi-vn)	107

IRISpec CA/CB/CC

REF 800-7702

Rx Only in the U.S.A.

INTENDED USE

IRISpec CA/CB/CC is an assayed QC material for monitoring of urine chemistry analytes and devices as listed in this Instruction for Use.

For *in vitro* diagnostic use only.

WARNING AND PRECAUTIONS

- All biological materials used in this control set are of non-human origin.
- Appropriate precautions in the collection, handling storage and disposal of specimens should be exercised accordingly.
- Wear appropriate personal protective equipment when handling.
- Avoid skin and eye contact. Refer to the Safety Data Sheet for first aid measures.
- Consult local or state authorities with regard to recommended method of disposal.
- Do not pipette.

GHS HAZARD CLASSIFICATION

IRISpec CA Urine Chemistry DANGER
Control

H360

P201

P280

P308+P313

May damage fertility or the unborn child.
Obtain special instructions before use.
Wear protective gloves, protective clothing and eye/face protection.
IF exposed or concerned: Get medical advice/attention.
Sodium Borate Decahydrate 1 - 2%

IRISpec CB Urine Chemistry DANGER
Control

H360

EUH208
P201

P280

P308+P313

May damage fertility or the unborn child.
May produce an allergic reaction.
Obtain special instructions before use.
Wear protective gloves, protective clothing and eye/face protection.
IF exposed or concerned: Get medical advice/attention.
Boric Acid 0.1 - < 0.3%
Subtilisin < 0.1%

IRISpec CC Urine Chemistry DANGER
Control

H360

May damage fertility or the unborn child.

P201	Obtain special instructions before use.
P280	Wear protective gloves, protective clothing and eye/face protection.
P308+P313	IF exposed or concerned: Get medical advice/attention. Boric Acid 0.3 - 1%

Safety Data Sheet is available at beckmancoulter.com/techdocs

REACTIVE INGREDIENTS

- Glucose, 11 g/L
- Bilirubin (bovine), 80 mg/L
- Sodium nitrite, 6.7 mg/L
- Ascorbic acid, 100 mg/dL
- Hemoglobin (bovine) 14.3 mg/L
- 2,5-dimethylindole, 54 mg/L
- Albumin (bovine) 4.5 g/L
- Calcium chloride, 8.0 g/L
- Acetone, 16 mL/L
- Microbial esterase, 120 mg/L

NON-REACTIVE INGREDIENTS

- Stabilizing Agents
- Buffers
- Salts
- Preservatives
- Deionized Water

STORAGE AND STABILITY

- Store at 2°C to 8°C
- Refrigerate immediately upon receipt
- IRISpec CA, CB and CC are stable for six months from the date of manufacture
- Open expiration: 15 days at 2°C to 8°C
- Protect from light
- Do not use solutions after the expiration date printed on the bottle labels

MATERIALS PROVIDED

- 3 x 100 mL IRISpec CA Control
- 3 x 100 mL IRISpec CB Control
- 3 x 100 mL IRISpec CC Control

MATERIALS REQUIRED, BUT NOT PROVIDED

- Urine Chemistry Strips
- Sample Tubes
- Personal Protective Equipment

PROCEDURE

1. Remove CA/CB/CC set from the refrigerator; pour an aliquot from each bottle into a separate tube. Ensure sufficient volume to wet the strip if testing manually or follow instrument manufacturer's instructions. Return tightly capped control bottles to the refrigerator immediately.
2. Allow aliquots to warm to room temperature (20° to 28°C). Use control aliquots within one hour of pouring. Protect aliquots from light.
3. Treat each aliquot as a separate control material. If using visual read strips, do not dip more than one strip into an aliquot.
4. Compare results obtained for each control with the Assay Value – Table of Acceptable Result Ranges (shown below) for the appropriate strip or strip-instrument combination. Ranges given are intended only as guidelines. More stringent ranges may be validated by each laboratory. However, establishing more lenient ranges is not acceptable. Note: a value less than 40mg/dL for ascorbic acid should prompt the operator to troubleshoot the chemistry strip, the control or instrument (if applicable).
5. Discard aliquots after single use. Do not pour back into bottle.

ASSAY VALUES – TABLE OF ACCEPTABLE RESULT RANGES			
iChemVELOCITY and DxU Urine Chemistry Systems with iChemVELOCITY Urine Chemistry Strips			
Analyte	IRISpec CA	IRISpec CB	IRISpec CC
Bilirubin	0.5 – ≥ 8.0 mg/dL 8.5 – ≥ 140 µmol/L (+) – (+++)	NEG.	NEG.
Urobilinogen	NORM.	2 – ≥ 8 mg/dL 34 – ≥ 140 µmol/L (+) – (+++)	NORM.
Ketone	NEG.	10 – > 100 mg/dL 1 – > 10 mmol/L (+) – (+++)	NEG.
Ascorbic Acid	NEG.	NEG.	40 mg/dL 2.28 mmol/L ++
Glucose	200 – ≥ 1000 mg/dL 11 – ≥ 56 mmol/L (++) – (++++)	NA	NEG.
Protein	100 – ≥ 600 mg/dL 1 – ≥ 6 g/L (++) – (++++)	NEG.	NEG.
Blood	0.03 – 0.5 mg/dL .3 – 5 mg/L Trace – (++)	NEG.	NEG.
pH	8 – 9	5 – 6	5 – 6
Nitrite	NEG.	0.1 – 0.2 mg/dL (+) – (++)	NEG.
Leukocytes	NEG.	75 – 500 WBC/µL (+) – (+++)	NEG.
SG*	1.005 – 1.020	1.010 – 1.020	NA

LIMITATION OF THE PROCEDURE

Minimize oxygen exposure of product. Do not shake bottles. Improper storage and handling may reduce analyte concentration (specifically Bilirubin) and may prematurely result in a negative reading. A color change from orange to green for the CA control may indicate such a process has occurred. The CB control may become turbid with time; this does not impact the performance. Each laboratory should establish its own precision parameters.

TECHNICAL SUPPORT

For addition information, or if damaged product is received, call Beckman Coulter Technical Support at 800-854-3633 (USA or Canada) or contact your local Beckman Coulter Representative.

A Glossary of Symbols (PN C27648) is available at beckmancoulter.com/techdocs.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

REVISION HISTORY

REVISION	DATE	DESCRIPTION
AA	December 2018	Initial Release

REVISION	DATE	DESCRIPTION
AB	August 2021	Updated Proprietary Name
		Updated website address
		Updated Legal Manufacturer Phone Number
		Updated copyright year.

REVISION	DATE	DESCRIPTION
AC	December 2023	Updated PROCEDURE section
		Updated Brazilian Legal Manufacturer Address
		Updated copyright year.

REVISION	DATE	DESCRIPTION
AD	February 2025	Updated Assay Table

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Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mills, Co. Clare, Ireland +(353) (0) 65 683 1100