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**Declaration of Conformity in accordance with the In Vitro
Diagnostic Directive 98/79/EC**

1. Manufacturer and or Authorised Representatives Details

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle upon Tyne
NE12 8EW
United Kingdom

Tel: 0191 215 0567
Fax: 0191 215 1152

2. The following Devices, classified as 'all other IVD Medical Devices', conform to the relevant provisions of the IVD Directive 98/79/EC and are CE marked in accordance with Annex III.

Bond™ System Reagents for use on the automated Bond system. As listed in
Device Schedule A

3. Description of the device

Bond™ Reagents for use on the automated Bond System.

4. Signed on behalf of the manufacturer by

Signed:  Date: 31 OCT 2018
Name: JEREMY ALLEN
Job Title: SNR QRA MANAGER

Device Schedule A (Continued)

Device Schedule A

The Bond™ - System Reagents includes the following products:

System Name and IFU Name	Catalogue Number
Bond™ Polymer Define Detection	DS9713
Bond™ Polymer Intense Detection	DS9958
Bond™ Polymer Refine Detection	DS9800
Bond™ Polymer Refine Red Detection	DS9390
ChromoPlex™ 1 Dual Detection	DS9477
ChromoPlex™ 1 Dual Detection – 50 Test	DS9665

Cleaning System Name and IFU Name	Catalogue Number
Bond™ Aspirating Probe Cleaning System	CS9100

Bond ISH Probes	Product Code	IFU Product Name
Bond™ Anti-Biotin, 7.5mL	PB0589	Bond™ Ready-to-Use ISH EBER Probe
Bond™ Anti-Fluorescein Antibody, 3.75 mL	PB0645	Bond™ Ready-to-Use ISH Kappa Probe
Bond™ Anti-Biotin, 7.5mL	PB0669	Bond™ Ready-to-Use ISH Lambda Probe
Bond™ Anti-Fluorescein Antibody, 3.75 mL	PB0785	Bond™ Ready-to-Use ISH RNA Positive Control Probe
Bond™ Anti-Biotin, 7.5mL	PB0809	Bond™ Ready-to-Use ISH RNA Negative Control Probe
Bond™ Anti-Fluorescein Antibody, 3.75 mL	PB0731	Bond™ Ready-To-Use ISH DNA Negative Control
Bond™ Anti-Biotin, 7.5mL	PB0682	Bond™ Ready-To-Use ISH DNA Positive Control
Bond™ Anti-Fluorescein Antibody, 3.75 mL	PB0780	Bond™ Ready-To-Use ISH HPV Probe (Subtypes 6, 11)
Bond™ Anti-Biotin, 7.5mL	PB0829	Bond™ Ready-to-Use ISH HPV Probe (Subtypes 16, 18, 31, 33, 51)

Bond Ancillary Reagents	Product Code	IFU Product Name
Bond™ Anti-Biotin, 7.5mL	AR0584	Bond™ Ready-To-Use Anti-Biotin Antibody
Bond™ Anti-Fluorescein Antibody, 3.75 mL	AR0833	Bond™ Ready-to-Use Anti-Fluorescein Antibody 3.75 mL
Bond™ Anti-Fluorescein Antibody 15mL	AR0222	Bond™ Ready-to-Use Anti-Fluorescein Antibody 15 mL
Bond™ DAB Enhancer	AR9432	Bond™ DAB Enhancer
Bond™ Dewax Solution	AR9222	Bond™ Dewax Solution
Bond™ Epitope Retrieval Solution 1	AR9961	Bond™ Epitope Retrieval Solution 1
Bond™ Epitope Retrieval Solution 2	AR9640	Bond™ Epitope Retrieval Solution 2

Device Schedule A (Continued)

Bond Ancillary Reagents	Product Code	IFU Product Name
Bond™ Enzyme Pretreatment Kit	AR9551	Bond™ Enzyme Pretreatment Kit
BOND Hybridization Solution	AR9037	BOND Hybridization Solution
Bond™ Primary Antibody Diluent, 500 mL	AR9352	Bond™ Primary Antibody Diluent
Bond™ Stringency Wash Solution, 3.75mL	AR0633	Bond™ Ready-To-Use Stringency Wash Solution
Bond™ Wash Solution 10X Concentrate, 1 L	AR9590	Bond™ Wash Solution 10X Concentrate

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ Alpha-Methylacyl-CoA Racemase (EPMU1), 7mL	PA0210	Bond™ Ready-to-Use Primary Antibody Alpha Methylacyl-CoA Racemase (EPMU1)
Bond™ Alpha-Fetoprotein (C3), 7 mL	PA0963	Bond™ Ready-To-Use Primary Antibody Alpha Fetoprotein (C3)
Bond™ Anaplastic Lymphoma Kinase (5A4), 7 mL	PA0306	Bond™ Ready-to-Use Primary Antibody Anaplastic Lymphoma Kinase (5A4)
Bond™ Anaplastic Lymphoma Kinase (5A4), 7mL	PA0831	Bond™ Ready-to-Use Primary Antibody Anaplastic Lymphoma Kinase (5A4)
Bond™ Bcl-2 (bcl-2/100/D5), 7 mL	PA0117	Bond™ Ready-to-Use Primary Antibody Bcl-2 (bcl-2/100/D5)
Bond™ Bcl-6 (LN22), 7 mL	PA0204	Bond™ Ready-to-Use Primary Antibody Bcl-6 (LN22)
Bond™ B cell Specific Octamer Binding Protein-1 (TG14), 7 mL	PA0558	Bond™ Ready-To-Use Primary Antibody B Cell Specific Octamer Binding Protein-1 (TG14)
Bond™ Beta-Catenin (17C2), 7mL	PA0083	Bond™ Ready-to-Use Primary Antibody Beta-Catenin (17C2)
Bond™ CA19-9 (Sialyl Lewis ^a) (C24:5:1:4), 7 mL	PA0424	Bond™ Ready-to-Use Primary Antibody CA19-9 (Sialyl Lewis ^a) (C241:5:1:4)
Bond™ CA125 (Ov185:1), 7 mL	PA0539	Bond™ Ready-to-Use Primary Antibody CA125 (Ovarian Cancer Antigen) (Ov185:1)
Bond™ Calcitonin (Polyclonal), 7 mL	PA0406	Bond™ Ready-To-Use Primary Antibody Calcitonin (Polyclonal)
Bond™ Calponin (Basic) (26A11), 7 mL	PA0416	Bond™ Ready-To-Use Primary Antibody Calponin (Basic) (26A11)
Bond™ Calretinin (CAL6), 7 mL	PA0346	Bond™ Ready-to-Use Primary Antibody Calretinin (CAL6)
Bond™ Carcinoembryonic Antigen (CD66e) (COL-1), 7 mL	PA0848	Bond™ Ready-to-Use Primary Antibody Carcinoembryonic Antigen (CD66e) (COL-1)
Bond™ Carcinoembryonic Antigen (II-7), 7 mL	PA0004	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Carcinoembryonic Antigen (II-7)
Bond™ c-erbB-2 Oncoprotein (CB11), 13.5mL	PA0571	Bond™ Ready-to-Use Primary Antibody c-erbB-2 Oncoprotein
Bond™ c-erbB-2 Oncoprotein (CB11), 7mL	PA0983	Bond™ Ready-to-Use Primary Antibody c-erbB-2 Oncoprotein
Bond™ CD1a (MTB1), 7 mL	PA0235	Bond™ Ready-to-Use Primary Antibody CD1a (MTB1)
Bond™ CD2 (11F11), 7 mL	PA0271	Bond™ Ready-to-Use Primary Antibody CD2 (11F11)
Bond™ CD3 (LN10), 7 mL	PA0553	Bond™ Ready-to-Use Primary Antibody CD3 (LN10)
Bond™ CD3 (LN10), 30 mL	PA0122	Bond™ Ready-to-Use Primary Antibody CD3 (LN10)
Bond™ CD4 (4B12), 7 mL	PA0368	Bond™ Ready-to-Use Primary Antibody CD4 (4B12)
Bond™ CD4 (4B12), 7 mL	PA0427	Bond™ Ready-to-Use Primary Antibody CD4 (4B12)
Bond™ CD5 (4C7), 7 mL	PA0168	Bond™ Ready-to-Use Primary Antibody CD5 (4C7)
Bond™ CD7 (LP15), 7 mL	PA0266	Bond™ Ready-to-Use Primary Antibody CD7 (LP15)
Bond™ CD8 (4B11), 7 mL	PA0183	Bond™ Ready-to-Use Primary Antibody CD8 (4B11)
Bond™ CD10 (56C6), 30 mL	PA0131	Bond™ Ready-to-Use Primary Antibody CD10 (56C6)
Bond™ CD10 (56C6), 7 mL	PA0270	Bond™ Ready-to-Use Primary Antibody CD10 (56C6)
Bond™ CD11c (5D11), 7 mL	PA0554	Bond™ Ready-to-Use Primary Antibody CD11c (5D11)
Bond™ CD13 (38C12), 7 mL	PA0304	Bond™ Ready-to-Use Primary Antibody CD13 (38C12)
Bond™ CD15 (Carb-1), 7 mL	PA0039	Bond™ Ready-To-Use Primary Antibody CD15 (Carb-1)
Bond™ CD15 (MMA), 7 mL	PA0473	Bond™ Ready-to-Use Primary Antibody CD15 (MMA)
Bond™ CD19 (BT51E), 7 mL	PA0843	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		CD19 (BT51E)
Bond™ CD20 (L26), 7 mL	PA0200	Bond™ Ready-to-Use Primary Antibody CD20 (L26)
Bond™ CD20 (L26), 30 mL	PA0359	Bond™ Ready-to-Use Primary Antibody CD20 (L26)
Bond™ CD20 (MJ1), 7 mL	PA0906	Bond™ Ready-to-Use Primary Antibody CD20 (MJ1)
Bond™ CD21 (2G9), 7 mL	PA0171	Bond™ Ready-to-Use Primary Antibody CD21 (2G9)
Bond™ CD22 (FPC1), 7 mL	PA0249	Bond™ Ready-to-Use Primary Antibody CD22 (FPC1)
Bond™ CD23 (1B12), 7 mL	PA0169	Bond™ Ready-to-Use Primary Antibody CD23 (1B12)
Bond™ CD25 (4C9), 7 mL	PA0305	Bond™ Ready-to-Use Primary Antibody CD25 (4C9)
Bond™ CD30 (1G12), 7 mL	PA0153	Bond™ Ready-to-Use Primary Antibody CD30 (1G12)
Bond™ CD30 (JCM182), 7 mL	PA0790	Bond™ Ready-to-Use Primary Antibody CD30 (JCM182)
Bond™ CD31 (1A10), 7 mL	PA0250	Bond™ Ready-to-Use Primary Antibody CD31 (1A10)
Bond™ CD31 (JC70A), 7mL	PA0414	Bond™ Ready-to-Use Primary Antibody CD31 (JC70A)
Bond™ CD33 (PWS44), 7 mL	PA0555	Bond™ Ready-to-Use Primary Antibody CD33 (PWS44)
Bond™ CD34 (QBEnd/10), 7 mL	PA0212	Bond™ Ready-to-Use Primary Antibody CD34 (QBEnd/10)
Bond™ CD34 (QBEnd/10), 30 mL	PA0354	Bond™ Ready-to-Use Primary Antibody CD34 (QBEnd/10)
Bond™ CD43 (MT1), 7 mL	PA0938	Bond™ Ready-to-Use Primary Antibody CD43 (MT1)
Bond™ CD45 (X16/99), 7 mL	PA0042	Bond™ Ready-to-Use Primary Antibody CD45 (X16/99)
Bond™ CD45RO (UCHL1), 7 mL	PA0146	Bond™ Ready-To-Use Primary Antibody CD45RO (UCHL1)
Bond™ CD56 (CD564), 7 mL	PA0191	Bond™ Ready-to-Use Primary Antibody CD56 (CD564)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ CD57 (NK-1), 7 mL	PA0443	Bond™ Ready-To-Use Primary Antibody CD57 (NK-1)
Bond™ CD61 (2f2), 7 mL	PA0308	Bond™ Ready-to-Use Primary Antibody CD61 (2f2)
Bond™ CD68 (514H12), 7 mL	PA0273	Bond™ Ready-to-Use Primary Antibody CD68 (514H12)
Bond™ CD79a (11E3), 7 mL	PA0192	Bond™ Ready-to-Use Primary Antibody CD79a (11E3)
Bond™ CD79a (JCB117), 7 mL	PA0599	Bond™ Ready-to-Use Primary Antibody CD79a (JCB117)
Bond™ CD99 (12E7), 7 mL	PA0509	Bond™ Ready-To-Use Primary Antibody CD99 (12E7)
Bond™ CD103 (EP206), 7 mL	PA0374	Bond™ Ready-to-Use Primary Antibody CD103 (EP206)
Bond™ CD117 (EP10), 7mL	PA0007	Bond™ Ready-to-Use Primary Antibody CD117 (EP10)
Bond™ CD138 (Syndecan-1) (MI15), 7 mL	PA0088	Bond™ Ready-to-Use Primary Antibody CD138 (MI15)
Bond™ CD163 (10D6), 7 mL	PA0090	Bond™ Ready-to-Use Primary Antibody CD163 (10D6)
Bond™ CDX2 (AMT28), 7 mL	PA0535	Bond™ Ready-to-Use Primary Antibody CDX2 (AMT28)
Bond™ CDX2 (EP25), 7 mL	PA0375	Bond™ Ready-to-Use Primary Antibody CDX2 (EP25)
Bond™ Chromogranin A (5H7), 7 mL	PA0430	Bond™ Ready-to-Use Primary Antibody Chromogranin A (5H7)
Bond™ Chromogranin A (5H7), 7 mL	PA0515	Bond™ Ready-to-Use Primary Antibody Chromogranin A (5H7)
Bond™ Cyclin D1 (EP12), 7 mL	PA0046	Bond™ Ready-to-Use Primary Antibody Cyclin D1 (EP12)
Bond™ Cytokeratin 5 (XM26), 7 mL	PA0468	Bond™ Ready-to-Use Primary Antibody Cytokeratin 5 (XM26)
Bond™ Cytokeratin 7 (RN7), 7 mL	PA0942	Bond™ Ready-to-Use Primary Antibody Cytokeratin 7 (RN7)
Bond™ Cytokeratin 7 (RN7), 30mL	PA0138	Bond™ Ready-to-Use Primary Antibody Cytokeratin 7 (RN7)
Bond™ Cytokeratin 8 (TS1), 7 mL	PA0567	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Cytokeratin 8 (TS1)
Bond™ Cytokeratin 8/18 (5D3), 7 mL	PA0067	Bond™ Ready-to-Use Primary Antibody Cytokeratin 8/18 (5D3)
Bond™ Cytokeratin 14 (LL002), 7mL	PA0074	Bond™ Ready-to-Use Primary Antibody Cytokeratin 14 (LL002)
Bond™ Cytokeratin 17 (E3), 7mL	PA0114	Bond™ Ready-to-Use Primary Antibody Cytokeratin 17 (E3)
Bond™ Cytokeratin 19 (b170), 7 mL	PA0799	Bond™ Ready-to-Use Primary Antibody Cytokeratin 19 (b170)
Bond™ Cytokeratin 20 (Ks20.8), 7 mL	PA0022	Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (Ks20.8)
Bond™ Cytokeratin 20 (Ks20.8), 30 mL	PA0037	Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (Ks20.8)
Bond™ Cytokeratin 20 (PW31), 7 mL	PA0918	Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (PW31)
Bond™ Cytokeratin (AE1/AE3), 7mL	PA0094	Bond™ Ready-to-Use Primary Antibody Cytokeratin (AE1/AE3)
Bond™ Cytokeratin HMW (34βE12), 7 mL	PA0134	Bond™ Ready-to-Use Primary Antibody Cytokeratin HMW (34βE12)
Bond™ Desmin (DE-R-11), 7 mL	PA0032	Bond™ Ready-to-Use Primary Antibody Desmin (DE-R-11)
Bond™ DOG-1 (K9), 7 mL	PA0219	Bond™ Ready-to-Use Primary Antibody DOG-1 (K9)
Bond™ E-Cadherin (36B5), 7 mL	PA0387	Bond™ Ready-To-Use Primary Antibody E-Cadherin (36B5)
Bond™ Enhancer of Zeste Homolog 2 (Drosophila) (EZH2) (6A10), 7mL	PA0575	Bond™ Ready-to-Use Primary Antibody Enhancer of Zeste Homolog 2 (Drosophila) (EZH2) (6A10)
Bond™ Epithelial Membrane Antigen (GP1.4), 7 mL	PA0035	Bond™ Ready-to-Use Primary Antibody Epithelial Membrane Antigen (GP1.4)
Bond™ Estrogen Receptor (6F11), 7 mL	PA0151	Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody for Bond™
Bond™ Estrogen Receptor (6F11), 30 mL	PA0009	Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody for Bond™
Bond™ Factor XIIIa (E890.1), 7 mL	PA0449	Bond™ Ready-To-Use Primary Antibody Factor XIIIa (E980.1)
Bond™ Fascin (IM20), 7 mL	PA0420	Bond™ Ready-to-Use Primary Antibody Fascin (IM20)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ Galectin-3 (9C4), 7mL	PA0238	Bond™ Ready-to-Use Primary Antibody Galectin-3 (9C4)
Bond™ Gastrin (Polyclonal), 7mL	PA0681	Bond™ Ready-to-Use Primary Antibody Gastrin (Polyclonal)
Bond™ Glial Fibrillary Acidic Protein (GA5), 7 mL	PA0026	Bond™ Ready-to-Use Primary Antibody Glial Fibrillary Acidic Protein (GA5)
Bond™ Glucagon (Polyclonal), 7mL	PA0594	Bond™ Ready-to-Use Primary Antibody Glucagon (Polyclonal)
Bond™ Granzyme B (11F1), 7 mL	PA0291	Bond™ Ready-to-Use Primary Antibody Granzyme B (11F1)
Bond™ Gross Cystic Disease Fluid Protein-15 (23A3), 7 mL	PA0350	Bond™ Ready-to-Use Primary Antibody Gross Cystic Disease Fluid Protein-15 (23A3)
Bond™ Gross Cystic Disease Fluid Protein-15 (23A3), 7 mL	PA0708	Bond™ Ready-to-Use Primary Antibody Gross Cystic Disease Fluid Protein-15 (23A3)
Bond™ Human Chorionic Gonadotrophin (Polyclonal), 7 mL	PA0014	Bond™ Ready-To-Use Primary Antibody Human Chorionic Gonadotrophin (Polyclonal)
Bond™ Human Follicle Stimulating Hormone (INN-hFSH-60), 7 mL	PA0693	Bond™ Ready-to-Use Primary Antibody Human Follicle Stimulating Hormone (INN-hFSH-60)
Bond™ Human Growth Hormone (Polyclonal), 7 mL	PA0704	Bond™ Ready-to-Use Primary Antibody Human Growth Hormone (Polyclonal)
Bond™ Human Herpesvirus (Type 8) (13B10), 7mL	PA0050	Bond™ Ready-to-Use Primary Antibody Human Herpesvirus (Type 8)
Bond™ Immunoglobulin D (DRN1C), 7mL	PA0061	Bond™ Ready-to-Use Primary Antibody Immunoglobulin D (DRN1C)
Bond™ Immunoglobulin G (RWP49), 7mL	PA0905	Bond™ Ready-to-Use Primary Antibody Immunoglobulin G (RWP49)
Bond™ Immunoglobulin M (8H6), 7mL	PA0278	Bond™ Ready-to-Use Primary Antibody Immunoglobulin M (8H6)
Bond™ Inhibin Alpha (R1), 7 mL	PA0488	Bond™ Ready-to-Use Primary Antibody Inhibin Alpha (R1)
Bond™ Insulin (2D11-H5), 7 mL	PA0620	Bond™ Ready-to-Use Primary Antibody Insulin (2D11-H5)
Bond™ Ki67 (MM1), 7 mL	PA0118	Bond™ Ready-to-Use Primary Antibody Ki67 (MM1)
Bond™ Ki67 (MM1), 30 mL	PA0410	Bond™ Ready-to-Use Primary Antibody Ki67 (MM1)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ Ki67 (K2), 7 mL	PA0230	Bond™ Ready-to-Use Primary Antibody Ki67 (K2)
Bond™ Kappa Light Chain (CH15), 7mL	PA0606	Bond™ Ready-to-Use Primary Antibody Kappa Light Chain (CH15)
Bond™ Lambda Light Chain (SHL53), 7mL	PA0570	Bond™ Ready-to-Use Primary Antibody Lambda Light Chain (SHL53)
Bond™ Luteinizing Hormone (C93), 7 mL	PA0655	Bond™ Ready-to-Use Primary Antibody Luteinizing Hormone (C93)
Bond™ Macrophage Marker (MAC387), 7 mL	PA0752	Bond™ Ready-to-Use Primary Antibody Macrophage Marker (MAC387)
Bond™ Mammaglobin (EP249), 7 mL	PA0378	Bond™ Ready-to-Use Primary Antibody Mammaglobin (EP249)
Bond™ Mast Cell Tryptase (10D11), 7 mL	PA0019	Bond™ Ready-To-Use Primary Antibody Mast Cell Tryptase (10D11)
Bond™ Melan A (A103), 7 mL	PA0233	Bond™ Ready-to-Use Primary Antibody Melan A (A103)
Bond™ Melan A (A103), 30 mL	PA0044	Bond™ Ready-to-Use Primary Antibody Melan A (A103)
Bond™ Melanoma Marker (HMB45), 7mL	PA0027	Bond™ Ready-to-Use Primary Antibody Melanoma Marker (HMB45)
Bond™ Melanoma Marker (HMB45), 30 mL	PA0625	Bond™ Ready-to-Use Primary Antibody Melanoma Marker (HMB45)
Bond™ Mesothelin (5B2), 7 mL	PA0373	Bond™ Ready-to-Use Primary Antibody Mesothelin (5B2)
Bond™ MLH1 (Mismatch Repair Protein), (ES05), 7 mL	PA0610	Bond™ Ready-to-Use Primary Antibody MLH1 (Mismatch Repair Protein) (ES05)
Bond™ MSH2 (Mismatch Repair Protein) (25D12),7 mL	PA0048	Bond™ Ready-to-Use Primary Antibody MSH2 (Mismatch Repair Protein) (25D12)
Bond™ Muc-1 Glycoprotein (Ma695), 7mL	PA0051	Bond™ Ready-to-Use Primary Antibody Muc-1 Glycoprotein (Ma695)
Bond™ Muc-2 Glycoprotein (Ccp58), 7mL	PA0155	Bond™ Ready-to-Use Primary Antibody Muc-2 (Ccp58)
Bond™ Muc-5AC Glycoprotein (CLH2), 7mL	PA0052	Bond™ Ready-to-Use Primary Antibody Muc-5AC Glycoprotein (CLH2)
Bond™ Muc-6 Glycoprotein (CLH5), 7mL	PA0053	Bond™ Ready-to-Use Primary Antibody Muc-6 Glycoprotein (CLH5)
Bond™ Multi-Cytokeratin (AE1/AE3), 7 mL	PA0909	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Multi-Cytokeratin (AE1/AE3)
Bond™ Multi-Cytokeratin (AE1/AE3), 7 mL	PA0012	Bond™ Ready-to-Use Primary Antibody Multi-Cytokeratin (AE1/AE3)
Bond™ Multiple Myeloma Oncogene 1 (MUM1), (EAU32), 7 mL	PA0129	Bond™ Ready-to-Use Primary Antibody Multiple Myeloma Oncogene 1 (MUM1) (EAU32)
Bond™ Muramidase (Lysozyme) (Polyclonal), 7 mL	PA0391	Bond™ Ready-to-Use Primary Antibody Muramidase (Lysozyme) (Polyclonal)
Bond™ Muscle Specific Actin (HHF35), 7 mL	PA0258	Bond™ Ready-to-Use Primary Antibody Muscle Specific Actin (HHF35)
Bond™ Myeloperoxidase (59A5), 7 mL	PA0491	Bond™ Ready-to-Use Primary Antibody Myeloperoxidase (59A5)
Bond™ Myf-4 (Rhabdomyosarcoma Marker) (LO26), 7 mL	PA0226	Bond™ Ready-to-Use Primary Antibody Myf-4 (Rhabdomyosarcoma Marker) (LO26)
Bond™ Myoglobin (MYO18), 7 mL	PA0727	Bond™ Ready-to-Use Primary Antibody Myoglobin (MYO18)
Bond™ Myosin Heavy Chain (Smooth muscle) (S131), 7 mL	PA0493	Bond™ Ready-to-Use Primary Antibody Myosin Heavy Chain (smooth muscle) (S131)
Bond™ Napsin A (IP64), 7mL	PA0064	Bond™ Ready-to-Use Primary Antibody Napsin A (IP64)
Bond™ Negative (Rabbit), 7 mL	PA0777	Bond™ Ready-to-Use Negative Control Negative (Rabbit)
Bond™ Negative (Mouse), 7 mL	PA0996	Bond™ Ready-to-Use Negative Control Negative (Mouse)
Bond™ Neurofilament 200kD (N52.1.7), 7 mL	PA0371	Bond™ Ready-to-Use Primary Antibody Neurofilament 200 kD (N52.1.7)
Bond™ Neuron Specific Enolase (22C9), 7 mL	PA0435	Bond™ Ready-to-Use Primary Antibody Neuron Specific Enolase (22C9)
Bond™ Oct-2 (Oct-207), 7 mL	PA0532	Bond™ Ready-to-Use Primary Antibody Oct-2 (Oct-207)
Bond™ Oct-3/4 (N1NK), 7 mL	PA0934	Bond™ Ready-to-Use Primary Antibody Oct-3/4 (N1NK)
Bond™ Oct-3/4 (N1NK), 7 mL	PA0193	Bond™ Ready-to-Use Primary Antibody Oct-3/4 (N1NK)
Bond™ p120 Catenin (EP66), 7 mL	PA0379	Bond™ Ready-to-Use Primary Antibody p120 Catenin (EP66)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ p53 (DO-7), 7 mL	PA0057	Bond™ Ready-to-Use Primary Antibody p53 (DO-7)
Bond™ p63 (7JUL), 7 mL	PA0103	Bond™ Ready-to-Use Primary Antibody p63 (7JUL)
Bond™ Pax-5 (1EW), 7 mL	PA0552	Bond™ Ready-To-Use Primary Antibody Pax-5 (1EW)
Bond™ Placental Alkaline Phosphatase (8A9), 7 mL	PA0161	Bond™ Ready-to-Use Primary Antibody Placental Alkaline Phosphatase (8A9)
Bond™ Progesterone Receptor (16), 7 mL	PA0312	Bond™ Ready-to-Use Primary Antibody Progesterone Receptor (16)
Bond™ Programmed Death Ligand 1 (73-10), 7mL	PA0832	Bond™ Ready-to-Use Primary Antibody Programmed Death Ligand 1 (73-10)
Bond™ Prostatic Acid Phosphatase (PASE/4LJ), 7 mL	PA0006	Bond™ Ready-to-Use Primary Antibody Prostatic Acid Phosphatase (PASE/4LJ)
Bond™ Protein Gene Product 9.5 (10A1), 7 mL	PA0286	Bond™ Ready-to-Use Primary Antibody Protein Gene Product 9.5 (10A1)
Bond™ S-100 (Polyclonal), 7 mL	PA0900	Bond™ Ready-to-Use Primary Antibody S-100 (Polyclonal)
Bond™ Serotonin (Polyclonal), 7 mL	PA0736	Bond™ Ready-to-Use Primary Antibody Serotonin (Polyclonal)
Bond™ Smooth Muscle Actin (asm-1), 7 mL	PA0943	Bond™ Ready-to-Use Primary Antibody Smooth Muscle Actin (alpha sm-1)
Bond™ Somatostatin (Polyclonal), 7 mL	PA0331	Bond™ Ready-To-Use Primary Antibody Somatostatin (Polyclonal)
Bond™ Synaptophysin (27G12), 7 mL	PA0299	Bond™ Ready-to-Use Primary Antibody Synaptophysin (27G12)
Bond™ Tartrate-Resistance Acid Phosphatase (26E5), 7ml	PA0093	Bond™ Ready-to-Use Primary Antibody Tartrate-Resistant Acid Phosphatase (26E5)
Bond™ Terminal Deoxynucleotidyl Transferase (SEN28), 7 mL	PA0339	Bond™ Ready-to-Use Primary Antibody Terminal Deoxynucleotidyl Transferase (SEN28)
Bond™ Thyroglobulin (1D4), 7 mL	PA0025	Bond™ Ready-To-Use Primary Antibody Thyroglobulin (1D4)
Bond™ Thyroid Stimulating Hormone (QB2/6),7 mL	PA0776	Bond™ Ready-to-Use Primary Antibody Thyroid Stimulating Hormone (QB2/6)
Bond™ Thyroid Transcription Factor-1 (SPT24), 7 mL	PA0364	Bond™ Ready-to-Use Primary Antibody Thyroid Transcription Factor-1 (SPT24)
Bond™ Tyrosinase (T311), 7 mL	PA0322	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Tyrosinase (T311)
Bond™ Villin (CWWB1), 7 mL	PA0106	Bond™ Ready-to-Use Primary Antibody Villin (CWWB1)
Bond™ Vimentin (SRL33), 7 mL	PA0033	Bond™ Ready-To-Use Primary Antibody Vimentin (SRL33)
Bond™ Vimentin (V9) 7mL	PA0640	Bond™ Ready-to-Use Primary Antibody Vimentin (V9)
Bond™ von Willebrand Factor (36B11), 7 mL	PA0400	Bond™ Ready-to-Use Primary Antibody von Willebrand Factor (36B11)
Bond™ von Willebrand Factor (36B11), 7 mL	PA0055	Bond™ Ready-to-Use Primary Antibody von Willebrand Factor (36B11)
Bond™ Wilms' Tumor (WT49), 7 mL	PA0562	Bond™ Ready-to-Use Primary Antibody Wilms' Tumor (WT49)
Bond™ ZAP-70 (L453R), 7 mL	PA0998	Bond™ Ready-to-Use Primary Antibody ZAP-70 (L453R)

Consumables	Product Code
Bond™ Open Containers 7 ml	OP79193
Bond™ Open Containers 30 ml	OP309700
Bond™ Titration Kit	OPT9049
Bond™ Titration Container Inserts – 50 Pack	OPT9719

Signed: 

Date: 31 OCT 2018

Name: JEREMY ALLEN

Job Title: SNR QRA MANAGER

**Declaration of Conformity In accordance with the In Vitro
Diagnostic Directive 98/79/EC**

1. Manufacturer and or Authorised Representatives Details

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle upon Tyne
NE12 8EW
United Kingdom

Tel: 0191 215 0567
Fax: 0191 215 1152

2. The following Devices, classified as 'all other IVD Medical Devices', conform to the relevant provisions of the IVD Directive 98/79/EC and are CE marked in accordance with Annex III.

Bond™ System Reagents for use on the automated Bond system. As listed in
Device Schedule A

3. Description of the device

Bond™ Reagents for use on the automated Bond System.

4. Signed on behalf of the manufacturer by

Signed: *P. Lloyd*

Date: *23 NOV 2016*

Name: *P. Lloyd*

Job Title: *QA Manager*

Device Schedule A (Continued)

Device Schedule A

The Bond™ - System Reagents includes the following products:

System Name and IFU Name	Catalogue Number
Bond™ Polymer Define Detection	DS9713
Bond™ Polymer Intense Detection	DS9958
Bond™ Polymer Refine Detection	DS9800
Bond™ Polymer Refine Red Detection	DS9390
ChromoPlex™ 1 Dual Detection	DS9477
ChromoPlex™ 1 Dual Detection – 50 Test	DS9665

Cleaning System Name and IFU Name	Catalogue Number
Bond™ Aspirating Probe Cleaning System	CS9100

Bond ISH Probes	Product Code	IFU Product Name
EBER Probe, 5.5 mL	PB0589	Bond™ Ready-to-Use ISH EBER Probe
Kappa Probe, 5.5 mL	PB0645	Bond™ Ready-to-Use ISH Kappa Probe
Lambda Probe, 5.5 mL	PB0669	Bond™ Ready-to-Use ISH Lambda Probe
RNA Positive Control Probe, 5.5 mL	PB0785	Bond™ Ready-to-Use ISH RNA Positive Control Probe
RNA Negative Control Probe, 5.5 mL	PB0809	Bond™ Ready-to-Use ISH RNA Negative Control Probe
DNA Negative Control 6.25 ml	PB0731	Bond™ Ready-To-Use ISH DNA Negative Control
DNA Positive Control 6.25 ml	PB0682	Bond™ Ready-To-Use ISH DNA Positive Control
HPV Probe (subtypes 6, 11) 6.25 ml	PB0780	Bond™ Ready-To-Use ISH HPV Probe (Subtypes 6, 11)
HPV Probe (subtypes 16, 18, 31, 33, 51) 6.25 ml	PB0829	Bond™ Ready-to-Use ISH HPV Probe (Subtypes 16, 18, 31, 33, 51)

Bond Ancillary Reagents	Product Code	IFU Product Name
Bond™ Anti-Biotin, 7.5mL	AR0584	Bond™ Ready-To-Use Anti-Biotin Antibody
Bond™ Anti-Fluorescein Antibody, 3.75 mL	AR0833	Bond™ Ready-to-Use Anti-Fluorescein Antibody 3.75 mL
Bond™ Anti-Fluorescein Antibody 15mL	AR0222	Bond™ Ready-to-Use Anti-Fluorescein Antibody 15 mL
Bond™ DAB Enhancer	AR9432	Bond™ DAB Enhancer
Bond™ Dewax Solution	AR9222	Bond™ Dewax Solution
Bond™ Epitope Retrieval Solution 1	AR9961	Bond™ Epitope Retrieval Solution 1
Bond™ Epitope Retrieval Solution 2	AR9640	Bond™ Epitope Retrieval Solution 2

Device Schedule A (Continued)

Bond Ancillary Reagents	Product Code	IFU Product Name
Bond™ Enzyme Pretreatment Kit	AR9551	Bond™ Enzyme Pretreatment Kit
BOND Hybridization Solution	AR9037	BOND Hybridization Solution
Bond™ Primary Antibody Diluent, 500 mL	AR9352	Bond™ Primary Antibody Diluent
Bond™ Stringency Wash Solution, 3.75mL	AR0633	Bond™ Ready-To-Use Stringency Wash Solution
Bond™ Wash Solution 10X Concentrate, 1 L	AR9590	Bond™ Wash Solution 10X Concentrate

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ Alpha-Methylacyl-CoA Racemase (EPMU1), 7mL	PA0210	Bond™ Ready-To-Use Primary Antibody Alpha Methylacyl-CoA Racemase (EPMU1)
Bond™ Alpha-Fetoprotein (C3), 7 mL	PA0963	Bond™ Ready-To-Use Primary Antibody Alpha Fetoprotein (C3)
Bond™ Anaplastic Lymphoma Kinase (5A4), 7 mL	PA0306	Bond™ Ready-to-Use Primary Antibody Anaplastic Lymphoma Kinase (5A4)
Bond™ Bcl-2 (bcl-2/100/D5), 7 mL	PA0117	Bond™ Ready-to-Use Primary Antibody Bcl-2 (bcl-2/100/D5)
Bond™ Bcl-6 (LN22), 7 mL	PA0204	Bond™ Ready-to-Use Primary Antibody Bcl-6 (LN22)
Bond™ B cell Specific Octamer Binding Protein-1 (TG14), 7 mL	PA0558	Bond™ Ready-To-Use Primary Antibody B Cell Specific Octamer Binding Protein-1 (TG14)
Bond™ Beta-Catenin (17C2), 7mL	PA0083	Bond™ Ready-to-Use Primary Antibody Beta-Catenin (17C2)
Bond™ CA19-9 (Sialyl Lewis ^a) (C24:5:1:4), 7 mL	PA0424	Bond™ Ready-to-Use Primary Antibody CA19-9 (Sialyl Lewis ^a) (C24:5:1:4)
Bond™ CA125 (Ov185:1), 7 mL	PA0539	Bond™ Ready-To-Use Primary Antibody CA125 (Ovarian Cancer Antigen) (Ov185:1)
Bond™ Calcitonin (Polyclonal), 7 mL	PA0406	Bond™ Ready-To-Use Primary Antibody Calcitonin (Polyclonal)
Bond™ Calponin (Basic) (26A11), 7 mL	PA0416	Bond™ Ready-To-Use Primary Antibody Calponin (Basic) (26A11)
Bond™ Calretinin (CAL6), 7 mL	PA0346	Bond™ Ready-to-Use Primary Antibody Calretinin (CAL6)
Bond™ Carcinoembryonic Antigen (CD66e) (COL-1), 7 mL	PA0848	Bond™ Ready-to-Use Primary Antibody Carcinoembryonic Antigen (CD66e) (COL-1)
Bond™ Carcinoembryonic Antigen (II-7), 7 mL	PA0004	Bond™ Ready-to-Use Primary Antibody Carcinoembryonic Antigen (II-7)
Bond™ CD1a (MTB1), 7 mL	PA0235	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		CD1a (MTB1)
Bond™ CD2 (11F11), 7 mL	PA0271	Bond™ Ready-to-Use Primary Antibody CD2 (11F11)
Bond™ CD3 (LN10), 7 mL	PA0553	Bond™ Ready-to-Use Primary Antibody CD3 (LN10)
Bond™ CD3 (LN10), 30 mL	PA0122	Bond™ Ready-to-Use Primary Antibody CD3 (LN10)
Bond™ CD4 (4B12), 7 mL	PA0368	Bond™ Ready-to-Use Primary Antibody CD4 (4B12)
Bond™ CD4 (4B12), 7 mL	PA0427	Bond™ Ready-to-Use Primary Antibody CD4 (4B12)
Bond™ CD5 (4C7), 7 mL	PA0168	Bond™ Ready-to-Use Primary Antibody CD5 (4C7)
Bond™ CD7 (LP15), 7 mL	PA0266	Bond™ Ready-to-Use Primary Antibody CD7 (LP15)
Bond™ CD8 (4B11), 7 mL	PA0183	Bond™ Ready-to-Use Primary Antibody CD8 (4B11)
Bond™ CD10 (56C6), 30 mL	PA0131	Bond™ Ready-to-Use Primary Antibody CD10 (56C6)
Bond™ CD10 (56C6), 7 mL	PA0270	Bond™ Ready-to-Use Primary Antibody CD10 (56C6)
Bond™ CD11c (5D11), 7 mL	PA0554	Bond™ Ready-to-Use Primary Antibody CD11c (5D11)
Bond™ CD13 (38C12), 7 mL	PA0304	Bond™ Ready-To-Use Primary Antibody CD13 (38C12)
Bond™ CD15 (Carb-1), 7 mL	PA0039	Bond™ Ready-To-Use Primary Antibody CD15 (Carb-1)
Bond™ CD15 (MMA), 7 mL	PA0473	Bond™ Ready-To-Use Primary Antibody CD15 (MMA)
Bond™ CD19 (BT51E), 7 mL	PA0843	Bond™ Ready-to-Use Primary Antibody CD19 (BT51E)
Bond™ CD20 (L26), 7 mL	PA0200	Bond™ Ready-to-Use Primary Antibody CD20 (L26)
Bond™ CD20 (L26), 30 mL	PA0359	Bond™ Ready-to-Use Primary Antibody CD20 (L26)
Bond™ CD20 (MJ1), 7 mL	PA0906	Bond™ Ready-to-Use Primary Antibody CD20 (MJ1)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ CD21 (2G9), 7 mL	PA0171	Bond™ Ready-to-Use Primary Antibody CD21 (2G9)
Bond™ CD22 (FPC1), 7 mL	PA0249	Bond™ Ready-to-Use Primary Antibody CD22 (FPC1)
Bond™ CD23 (1B12), 7 mL	PA0169	Bond™ Ready-to-Use Primary Antibody CD23 (1B12)
Bond™ CD25 (4C9), 7 mL	PA0305	Bond™ Ready-to-Use Primary Antibody CD25 (4C9)
Bond™ CD30 (1G12), 7 mL	PA0153	Bond™ Ready-to-Use Primary Antibody CD30 (1G12)
Bond™ CD30 (JCM182), 7 mL	PA0790	Bond™ Ready-to-Use Primary Antibody CD30 (JCM182)
Bond™ CD31 (1A10), 7 mL	PA0250	Bond™ Ready-to-Use Primary Antibody CD31 (1A10)
Bond™ CD31 (JC70A), 7mL	PA0414	Bond™ Ready-to-Use Primary Antibody CD31 (JC70A)
Bond™ CD33 (PWS44), 7 mL	PA0555	Bond™ Ready-to-Use Primary Antibody CD33 (PWS44)
Bond™ CD34 (QBEnd/10), 7 mL	PA0212	Bond™ Ready-to-Use Primary Antibody CD34 (QBEnd/10)
Bond™ CD34 (QBEnd/10), 30 mL	PA0354	Bond™ Ready-to-Use Primary Antibody CD34 (QBEnd/10)
Bond™ CD43 (MT1), 7 mL	PA0938	Bond™ Ready-to-Use Primary Antibody CD43 (MT1)
Bond™ CD45 (X16/99), 7 mL	PA0042	Bond™ Ready-to-Use Primary Antibody CD45 (X16/99)
Bond™ CD45RO (UCHL1), 7 mL	PA0146	Bond™ Ready-To-Use Primary Antibody CD45RO (UCHL1)
Bond™ CD56 (CD564), 7 mL	PA0191	Bond™ Ready-to-Use Primary Antibody CD56 (CD564)
Bond™ CD57 (NK-1), 7 mL	PA0443	Bond™ Ready-To-Use Primary Antibody CD57 (NK-1)
Bond™ CD61 (2f2), 7 mL	PA0308	Bond™ Ready-to-Use Primary Antibody CD61 (2f2)
Bond™ CD68 (514H12), 7 mL	PA0273	Bond™ Ready-to-Use Primary Antibody CD68 (514H12)
Bond™ CD79a (11E3), 7 mL	PA0192	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		CD79a (11E3)
Bond™ CD79a (JCB117), 7 mL	PA0599	Bond™ Ready-to-Use Primary Antibody CD79a (JCB117)
Bond™ CD99 (12E7), 7 mL	PA0509	Bond™ Ready-To-Use Primary Antibody CD99 (12E7)
Bond™ CD103 (EP206), 7 mL	PA0374	Bond™ Ready-To-Use Primary Antibody CD103 (EP206)
Bond™ CD117 (EP10), 7mL	PA0007	Bond™ Ready-to-Use Primary Antibody CD117 (EP10)
Bond™ CD138 (Syndecan-1) (MI15), 7 mL	PA0088	Bond™ Ready-to-Use Primary Antibody CD138 (MI15)
Bond™ CD163 (10D6), 7 mL	PA0090	Bond™ Ready-to-Use Primary Antibody CD163 (10D6)
Bond™ CDX2 (AMT28), 7 mL	PA0535	Bond™ Ready-to-Use Primary Antibody CDX2 (AMT28)
Bond™ CDX2 (EP25), 7 mL	PA0375	Bond™ Ready-to-Use Primary Antibody CDX2 (EP25)
Bond™ Chromogranin A (5H7), 7 mL	PA0430	Bond™ Ready-to-Use Primary Antibody Chromogranin A (5H7)
Bond™ Chromogranin A (5H7), 7 mL	PA0515	Bond™ Ready-to-Use Primary Antibody Chromogranin A (5H7)
Bond™ Cyclin D1 (EP12), 7 mL	PA0046	Bond™ Ready-to-Use Primary Antibody Cyclin D1 (EP12)
Bond™ Cytokeratin 5 (XM26), 7 mL	PA0468	Bond™ Ready-to-Use Primary Antibody Cytokeratin 5 (XM26)
Bond™ Cytokeratin 7 (RN7), 7 mL	PA0942	Bond™ Ready-to-Use Primary Antibody Cytokeratin 7 (RN7)
Bond™ Cytokeratin 7 (RN7), 30mL	PA0138	Bond™ Ready-to-Use Primary Antibody Cytokeratin 7 (RN7)
Bond™ Cytokeratin 8 (TS1), 7 mL	PA0567	Bond™ Ready-to-Use Primary Antibody Cytokeratin 8 (TS1)
Bond™ Cytokeratin 8/18 (5D3), 7 mL	PA0067	Bond™ Ready-to-Use Primary Antibody Cytokeratin 8/18 (5D3)
Bond™ Cytokeratin 14 (LL002), 7mL	PA0074	Bond™ Ready-to-Use Primary Antibody Cytokeratin 14 (LL002)
Bond™ Cytokeratin 17 (E3), 7mL	PA0114	Bond™ Ready-to-Use Primary Antibody Cytokeratin 17 (E3)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ Cytokeratin 19 (b170), 7 mL	PA0799	Bond™ Ready-to-Use Primary Antibody Cytokeratin 19 (b170)
Bond™ Cytokeratin 20 (Ks20.8), 7 mL	PA0022	Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (Ks20.8)
Bond™ Cytokeratin 20 (Ks20.8), 30 mL	PA0037	Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (Ks20.8)
Bond™ Cytokeratin 20 (PW31), 7 mL	PA0918	Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (PW31)
Bond™ Cytokeratin HMW (34βE12), 7 mL	PA0134	Bond™ Ready-To-Use Primary Antibody Cytokeratin HMW (34βE12)
Bond™ Desmin (DE-R-11), 7 mL	PA0032	Bond™ Ready-to-Use Primary Antibody Desmin (DE-R-11)
Bond™ DOG-1 (K9), 7 mL	PA0219	Bond™ Ready-to-Use Primary Antibody DOG-1 (K9)
Bond™ E-Cadherin (36B5), 7 mL	PA0387	Bond™ Ready-To-Use Primary Antibody E-Cadherin (36B5)
Bond™ Enhancer of Zeste Homolog 2 (Drosophila) (EZH2) (6A10), 7mL	PA0575	Bond™ Ready-to-Use Primary Antibody Enhancer of Zeste Homolog 2 (Drosophila) (EZH2) (6A10)
Bond™ Epithelial Membrane Antigen (GP1.4), 7 mL	PA0035	Bond™ Ready-to-Use Primary Antibody Epithelial Membrane Antigen (GP1.4)
Bond™ Estrogen Receptor (6F11), 7 mL	PA0151	Estrogen Receptor Clone 6F11Ready-To-Use Primary Antibody for Bond™
Bond™ Estrogen Receptor (6F11), 30 mL	PA0009	Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody for Bond™
Bond™ Factor XIIIa (E890.1), 7 mL	PA0449	Bond™ Ready-To-Use Primary Antibody Factor XIIIa (E980.1)
Bond™ Fascin (IM20), 7 mL	PA0420	Bond™ Ready-to-Use Primary Antibody Fascin (IM20)
Bond™ Galectin-3 (9C4), 7mL	PA0238	Bond™ Ready-to-Use Primary Antibody Galectin-3 (9C4)
Bond™ Gastrin (Polyclonal), 7mL	PA0681	Bond™ Ready-to-Use Primary Antibody Gastrin (Polyclonal)
Bond™ Glial Fibrillary Acidic Protein (GA5), 7 mL	PA0026	Bond™ Ready-to-Use Primary Antibody Glial Fibrillary Acidic Protein (GA5)
Bond™ Glucagon (Polyclonal), 7mL	PA0594	Bond™ Ready-to-Use Primary Antibody Glucagon (Polyclonal)
Bond™ Granzyme B (11F1), 7 mL	PA0291	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Granzyme B (11F1)
Bond™ Gross Cystic Disease Fluid Protein-15 (23A3), 7 mL	PA0350	Bond™ Ready-to-Use Primary Antibody Gross Cystic Disease Fluid Protein-15 (23A3)
Bond™ Gross Cystic Disease Fluid Protein-15 (23A3), 7 mL	PA0708	Bond™ Ready-to-Use Primary Antibody Gross Cystic Disease Fluid Protein-15 (23A3)
Bond™ Human Chorionic Gonadotrophin (Polyclonal), 7 mL	PA0014	Bond™ Ready-To-Use Primary Antibody Human Chorionic Gonadotrophin (Polyclonal)
Bond™ Human Follicle Stimulating Hormone (INN-hFSH-60), 7 mL	PA0693	Bond™ Ready-to-Use Primary Antibody Human Follicle Stimulating Hormone (INN-hFSH-60)
Bond™ Human Growth Hormone (Polyclonal), 7 mL	PA0704	Bond™ Ready-to-Use Primary Antibody Human Growth Hormone (Polyclonal)
Bond™ Immunoglobulin D (DRN1C), 7mL	PA0061	Bond™ Ready-to-Use Primary Antibody Immunoglobulin D (DRN1C)
Bond™ Immunoglobulin G (RWP49), 7mL	PA0905	Bond™ Ready-to-Use Primary Antibody Immunoglobulin G (RWP49)
Bond™ Immunoglobulin M (8H6), 7mL	PA0278	Bond™ Ready-to-Use Primary Antibody Immunoglobulin M (8H6)
Bond™ Inhibin Alpha (R1), 7 mL	PA0488	Bond™ Ready-to-Use Primary Antibody Inhibin Alpha (R1)
Bond™ Insulin (2D11-H5), 7 mL	PA0620	Bond™ Ready-to-Use Primary Antibody Insulin (2D11-H5)
Bond™ Ki67 (MM1), 7 mL	PA0118	Bond™ Ready-to-Use Primary Antibody Ki67 (MM1)
Bond™ Ki67 (MM1), 30 mL	PA0410	Bond™ Ready-to-Use Primary Antibody Ki67 (MM1)
Bond™ Ki67 (K2), 7 mL	PA0230	Bond™ Ready-to-Use Primary Antibody Ki67 (K2)
Bond™ Kappa Light Chain (CH15), 7mL	PA0606	Bond™ Ready-to-Use Primary Antibody Kappa Light Chain (CH15)
Bond™ Lambda Light Chain (SHL53), 7mL	PA0570	Bond™ Ready-to-Use Primary Antibody Lambda Light Chain (SHL53)
Bond™ Luteinizing Hormone (C93), 7 mL	PA0655	Bond™ Ready-to-Use Primary Antibody Luteinizing Hormone (C93)
Bond™ Macrophage Marker (MAC387), 7 mL	PA0752	Bond™ Ready-to-Use Primary Antibody Macrophage Marker (MAC387)
Bond™ Mammaglobin (EP249), 7 mL	PA0378	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Mammaglobin (EP249)
Bond™ Mast Cell Tryptase (10D11), 7 mL	PA0019	Bond™ Ready-To-Use Primary Antibody Mast Cell Tryptase (10D11)
Bond™ Melan A (A103), 7 mL	PA0233	Bond™ Ready-to-Use Primary Antibody Melan A (A103)
Bond™ Melan A (A103), 30 mL	PA0044	Bond™ Ready-to-Use Primary Antibody Melan A (A103)
Bond™ Melanoma Marker (HMB45), 7mL	PA0027	Bond™ Ready-to-Use Primary Antibody Melanoma Marker (HMB45)
Bond™ Melanoma Marker (HMB45), 30 mL	PA0625	Bond™ Ready-to-Use Primary Antibody Melanoma Marker (HMB45)
Bond™ Mesothelin (5B2), 7 mL	PA0373	Bond™ Ready-to-Use Primary Antibody Mesothelin (5B2)
Bond™ MLH1 (Mismatch Repair Protein), (ES05), 7 mL	PA0610	Bond™ Ready-to-Use Primary Antibody MLH1 (Mismatch Repair Protein) (ES05)
Bond™ MSH2 (Mismatch Repair Protein) (25D12), 7 mL	PA0048	Bond™ Ready-to-Use Primary Antibody MSH2 (Mismatch Repair Protein) (25D12)
Bond™ Muc-2 Glycoprotein (Ccp58), 7mL	PA0155	Bond™ Ready-to-Use Primary Antibody Muc-2 Glycoprotein (Ccp58)
Bond™ Multi-Cytokeratin (AE1/AE3), 7 mL	PA0909	Bond™ Ready-to-Use Primary Antibody Multi-Cytokeratin (AE1/AE3)
Bond™ Multi-Cytokeratin (AE1/AE3), 7 mL	PA0012	Bond™ Ready-to-Use Primary Antibody Multi-Cytokeratin (AE1/AE3)
Bond™ Multiple Myeloma Oncogene 1 (MUM1), (EAU32), 7 mL	PA0129	Bond™ Ready-to-Use Primary Antibody Multiple Myeloma Oncogene 1 (MUM1) (EAU32)
Bond™ Muramidase (Lysozyme) (Polyclonal), 7 mL	PA0391	Bond™ Ready-to-Use Primary Antibody Muramidase (Lysozyme) (Polyclonal)
Bond™ Muscle Specific Actin (HHF35), 7 mL	PA0258	Bond™ Ready-to-Use Primary Antibody Muscle Specific Actin (HHF35)
Bond™ Myeloperoxidase (59A5), 7 mL	PA0491	Bond™ Ready-to-Use Primary Antibody Myeloperoxidase (59A5)
Bond™ Myf-4 (Rhabdomyosarcoma Marker) (LO26), 7 mL	PA0226	Bond™ Ready-to-Use Primary Antibody Myf-4 (Rhabdomyosarcoma Marker) (LO26)
Bond™ Myoglobin (MYO18), 7 mL	PA0727	Bond™ Ready-to-Use Primary Antibody Myoglobin (MYO18)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ Myosin Heavy Chain (Smooth muscle) (S131), 7 mL	PA0493	Bond™ Ready-to-Use Primary Antibody Myosin Heavy Chain (smooth muscle) (S131)
Bond™ Napsin A (IP64), 7mL	PA0064	Bond™ Ready-to-Use Napsin A (IP64)
Bond™ Negative (Rabbit), 7 mL	PA0777	Bond™ Ready-to-Use Negative Control Negative (Rabbit)
Bond™ Negative (Mouse), 7 mL	PA0996	Bond™ Ready-to-Use Negative Control Negative (Mouse)
Bond™ Neurofilament 200kD (N52.1.7), 7 mL	PA0371	Bond™ Ready-to-Use Primary Antibody Neurofilament 200 kD (N52.1.7)
Bond™ Neuron Specific Enolase (22C9), 7 mL	PA0435	Bond™ Ready-to-Use Primary Antibody Neuron Specific Enolase (22C9)
Bond™ Oct-2 (Oct-207), 7 mL	PA0532	Bond™ Ready-to-Use Primary Antibody Oct-2 (Oct-207)
Bond™ Oct-3/4 (N1NK), 7 mL	PA0934	Bond™ Ready-to-Use Primary Antibody Oct-3/4 (N1NK)
Bond™ Oct-3/4 (N1NK), 7 mL	PA0193	Bond™ Ready-to-Use Primary Antibody Oct-3/4 (N1NK)
Bond™ p120 Catenin (EP66), 7 mL	PA0379	Bond™ Ready-to-Use Primary Antibody p120 Catenin (EP66)
Bond™ p53 (DO-7), 7 mL	PA0057	Bond™ Ready-to-Use Primary Antibody p53 (DO-7)
Bond™ p63 (7JUL), 7 mL	PA0103	Bond™ Ready-to-Use Primary Antibody p63 (7JUL)
Bond™ Pax-5 (1EW), 7 mL	PA0552	Bond™ Ready-To-Use Primary Antibody Pax-5 (1EW)
Bond™ Placental Alkaline Phosphatase (8A9), 7 mL	PA0161	Bond™ Ready-to-Use Primary Antibody Placental Alkaline Phosphatase (8A9)
Bond™ Progesterone Receptor (16), 7 mL	PA0312	Bond™ Ready-to-Use Primary Antibody Progesterone Receptor (16)
Bond™ Prostatic Acid Phosphatase (PASE/4LJ), 7 mL	PA0006	Bond™ Ready-to-Use Primary Antibody Prostatic Acid Phosphatase (PASE/4LJ)
Bond™ Protein Gene Product 9.5 (10A1), 7 mL	PA0286	Bond™ Ready-to-Use Primary Antibody Protein Gene Product 9.5 (10A1)
Bond™ S-100 (Polyclonal), 7 mL	PA0900	Bond™ Ready-to-Use Primary Antibody S-100 (Polyclonal)
Bond™ Serotonin (Polyclonal), 7 mL	PA0736	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Serotonin (Polyclonal)
Bond™ Smooth Muscle Actin (asm-1), 7 mL	PA0943	Bond™ Ready-to-Use Primary Antibody Smooth Muscle Actin (asm-1)
Bond™ Somatostatin (Polyclonal), 7 mL	PA0331	Bond™ Ready-To-Use Primary Antibody Somatostatin (Polyclonal)
Bond™ Synaptophysin (27G12), 7 mL	PA0299	Bond™ Ready-to-Use Primary Antibody Synaptophysin (27G12)
Bond™ Tartrate-Resistance Acid Phosphatase (26E5), 7ml	PA0093	Bond™ Ready-to-Use Primary Antibody Tartrate-Resistant Acid Phosphatase (26E5)
Bond™ Terminal Deoxynucleotidyl Transferase (SEN28), 7 mL	PA0339	Bond™ Ready-to-Use Primary Antibody Terminal Deoxynucleotidyl Transferase (SEN28)
Bond™ Thyroglobulin (1D4), 7 mL	PA0025	Bond™ Ready-To-Use Primary Antibody Thyroglobulin (1D4)
Bond™ Thyroid Stimulating Hormone (QB2/6), 7 mL	PA0776	Bond™ Ready-to-Use Primary Antibody Thyroid Stimulating Hormone (QB2/6)
Bond™ Thyroid Transcription Factor-1 (SPT24), 7 mL	PA0364	Bond™ Ready-to-Use Primary Antibody Thyroid Transcription Factor-1 (SPT24)
Bond™ Tyrosinase (T311), 7 mL	PA0322	Bond™ Ready-to-Use Primary Antibody Tyrosinase (T311)
Bond™ Villin (CWWB1), 7 mL	PA0106	Bond™ Ready-to-Use Primary Antibody Villin (CWWB1)
Bond™ Vimentin (SRL33), 7 mL	PA0033	Bond™ Ready-To-Use Primary Antibody Vimentin (SRL33)
Bond™ Vimentin (V9) 7mL	PA0640	Bond™ Ready-to-Use Primary Antibody Vimentin (V9)
Bond™ von Willebrand Factor (36B11), 7 mL	PA0400	Bond™ Ready-to-Use Primary Antibody von Willebrand Factor (36B11)
Bond™ von Willebrand Factor (36B11), 7 mL	PA0055	Bond™ Ready-to-Use Primary Antibody von Willebrand Factor (36B11)
Bond™ Wilms' Tumor (WT49), 7 mL	PA0562	Bond™ Ready-to-Use Primary Antibody Wilms' Tumor (WT49)
Bond™ ZAP-70 (L453R), 7 mL	PA0998	Bond™ Ready-to-Use Primary Antibody ZAP-70 (L453R)

Device Schedule A (Continued)

Consumables	Product Code
Bond™ Open Containers 7 ml	OP79193
Bond™ Open Containers 30 ml	OP309700
Bond™ Titration Kit	OPT9049
Bond™ Titration Container Inserts – 50 Pack	OPT9719

Signed:..... *P. Mend*

Date:..... *23 NOV 2016*

Name: *P. Mend*

Job Title: *QA Manager*

**Declaration of Conformity in accordance with the In Vitro
Diagnostic Directive 98/79/EC**

1. Manufacturer and or Authorised Representatives Details

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle upon Tyne
NE12 8EW
United Kingdom

Tel: 0191 215 0567
Fax: 0191 215 1152

2. The following Devices, classified as 'all other IVD Medical Devices', conform to the relevant provisions of the IVD Directive 98/79/EC and are CE marked in accordance with Annex III.

Bond™ System Reagents for use on the automated Bond system. As listed in
Device Schedule A

3. Description of the device

Bond™ Reagents for use on the automated Bond System.

4. Signed on behalf of the manufacturer by

Signed: 

Name: JEREMY ALLEN

Job Title: SNR QRA MANAGER

Date: 31 OCT 2018

Device Schedule A (Continued)

Device Schedule A

The Bond™ - System Reagents includes the following products:

System Name and IFU Name	Catalogue Number
Bond™ Polymer Define Detection	DS9713
Bond™ Polymer Intense Detection	DS9958
Bond™ Polymer Refine Detection	DS9800
Bond™ Polymer Refine Red Detection	DS9390
ChromoPlex™ 1 Dual Detection	DS9477
ChromoPlex™ 1 Dual Detection – 50 Test	DS9665

Cleaning System Name and IFU Name	Catalogue Number
Bond™ Aspirating Probe Cleaning System	CS9100

Bond ISH Probes	Product Code	IFU Product Name
Bond™ Anti-Biotin, 7.5mL	PB0589	Bond™ Ready-to-Use ISH EBER Probe
Bond™ Anti-Fluorescein Antibody, 3.75 mL	PB0645	Bond™ Ready-to-Use ISH Kappa Probe
Bond™ Anti-Biotin, 7.5mL	PB0669	Bond™ Ready-to-Use ISH Lambda Probe
Bond™ Anti-Fluorescein Antibody, 3.75 mL	PB0785	Bond™ Ready-to-Use ISH RNA Positive Control Probe
Bond™ Anti-Biotin, 7.5mL	PB0809	Bond™ Ready-to-Use ISH RNA Negative Control Probe
Bond™ Anti-Fluorescein Antibody, 3.75 mL	PB0731	Bond™ Ready-To-Use ISH DNA Negative Control
Bond™ Anti-Biotin, 7.5mL	PB0682	Bond™ Ready-To-Use ISH DNA Positive Control
Bond™ Anti-Fluorescein Antibody, 3.75 mL	PB0780	Bond™ Ready-To-Use ISH HPV Probe (Subtypes 6, 11)
Bond™ Anti-Biotin, 7.5mL	PB0829	Bond™ Ready-to-Use ISH HPV Probe (Subtypes 16, 18, 31, 33, 51)

Bond Ancillary Reagents	Product Code	IFU Product Name
Bond™ Anti-Biotin, 7.5mL	AR0584	Bond™ Ready-To-Use Anti-Biotin Antibody
Bond™ Anti-Fluorescein Antibody, 3.75 mL	AR0833	Bond™ Ready-to-Use Anti-Fluorescein Antibody 3.75 mL
Bond™ Anti-Fluorescein Antibody 15mL	AR0222	Bond™ Ready-to-Use Anti-Fluorescein Antibody 15 mL
Bond™ DAB Enhancer	AR9432	Bond™ DAB Enhancer
Bond™ Dewax Solution	AR9222	Bond™ Dewax Solution
Bond™ Epitope Retrieval Solution 1	AR9961	Bond™ Epitope Retrieval Solution 1
Bond™ Epitope Retrieval Solution 2	AR9640	Bond™ Epitope Retrieval Solution 2

Device Schedule A (Continued)

Bond Ancillary Reagents	Product Code	IFU Product Name
Bond™ Enzyme Pretreatment Kit	AR9551	Bond™ Enzyme Pretreatment Kit
BOND Hybridization Solution	AR9037	BOND Hybridization Solution
Bond™ Primary Antibody Diluent, 500 mL	AR9352	Bond™ Primary Antibody Diluent
Bond™ Stringency Wash Solution, 3.75mL	AR0633	Bond™ Ready-To-Use Stringency Wash Solution
Bond™ Wash Solution 10X Concentrate, 1 L	AR9590	Bond™ Wash Solution 10X Concentrate

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ Alpha-Methylacyl-CoA Racemase (EPMU1), 7mL	PA0210	Bond™ Ready-to-Use Primary Antibody Alpha Methylacyl-CoA Racemase (EPMU1)
Bond™ Alpha-Fetoprotein (C3), 7 mL	PA0963	Bond™ Ready-To-Use Primary Antibody Alpha Fetoprotein (C3)
Bond™ Anaplastic Lymphoma Kinase (5A4), 7 mL	PA0306	Bond™ Ready-to-Use Primary Antibody Anaplastic Lymphoma Kinase (5A4)
Bond™ Anaplastic Lymphoma Kinase (5A4), 7mL	PA0831	Bond™ Ready-to-Use Primary Antibody Anaplastic Lymphoma Kinase (5A4)
Bond™ Bcl-2 (bcl-2/100/D5), 7 mL	PA0117	Bond™ Ready-to-Use Primary Antibody Bcl-2 (bcl-2/100/D5)
Bond™ Bcl-6 (LN22), 7 mL	PA0204	Bond™ Ready-to-Use Primary Antibody Bcl-6 (LN22)
Bond™ B cell Specific Octamer Binding Protein-1 (TG14), 7 mL	PA0558	Bond™ Ready-To-Use Primary Antibody B Cell Specific Octamer Binding Protein-1 (TG14)
Bond™ Beta-Catenin (17C2), 7mL	PA0083	Bond™ Ready-to-Use Primary Antibody Beta-Catenin (17C2)
Bond™ CA19-9 (Sialyl Lewis ^a) (C24:5:1:4), 7 mL	PA0424	Bond™ Ready-to-Use Primary Antibody CA19-9 (Sialyl Lewis ^a) (C241:5:1:4)
Bond™ CA125 (Ov185:1), 7 mL	PA0539	Bond™ Ready-to-Use Primary Antibody CA125 (Ovarian Cancer Antigen) (Ov185:1)
Bond™ Calcitonin (Polyclonal), 7 mL	PA0406	Bond™ Ready-To-Use Primary Antibody Calcitonin (Polyclonal)
Bond™ Calponin (Basic) (26A11), 7 mL	PA0416	Bond™ Ready-To-Use Primary Antibody Calponin (Basic) (26A11)
Bond™ Calretinin (CAL6), 7 mL	PA0346	Bond™ Ready-to-Use Primary Antibody Calretinin (CAL6)
Bond™ Carcinoembryonic Antigen (CD66e) (COL-1), 7 mL	PA0848	Bond™ Ready-to-Use Primary Antibody Carcinoembryonic Antigen (CD66e) (COL-1)
Bond™ Carcinoembryonic Antigen (II-7), 7 mL	PA0004	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Carcinoembryonic Antigen (II-7)
Bond™ c-erbB-2 Oncoprotein (CB11), 13.5mL	PA0571	Bond™ Ready-to-Use Primary Antibody c-erbB-2 Oncoprotein
Bond™ c-erbB-2 Oncoprotein (CB11), 7mL	PA0983	Bond™ Ready-to-Use Primary Antibody c-erbB-2 Oncoprotein
Bond™ CD1a (MTB1), 7 mL	PA0235	Bond™ Ready-to-Use Primary Antibody CD1a (MTB1)
Bond™ CD2 (11F11), 7 mL	PA0271	Bond™ Ready-to-Use Primary Antibody CD2 (11F11)
Bond™ CD3 (LN10), 7 mL	PA0553	Bond™ Ready-to-Use Primary Antibody CD3 (LN10)
Bond™ CD3 (LN10), 30 mL	PA0122	Bond™ Ready-to-Use Primary Antibody CD3 (LN10)
Bond™ CD4 (4B12), 7 mL	PA0368	Bond™ Ready-to-Use Primary Antibody CD4 (4B12)
Bond™ CD4 (4B12), 7 mL	PA0427	Bond™ Ready-to-Use Primary Antibody CD4 (4B12)
Bond™ CD5 (4C7), 7 mL	PA0168	Bond™ Ready-to-Use Primary Antibody CD5 (4C7)
Bond™ CD7 (LP15), 7 mL	PA0266	Bond™ Ready-to-Use Primary Antibody CD7 (LP15)
Bond™ CD8 (4B11), 7 mL	PA0183	Bond™ Ready-to-Use Primary Antibody CD8 (4B11)
Bond™ CD10 (56C6), 30 mL	PA0131	Bond™ Ready-to-Use Primary Antibody CD10 (56C6)
Bond™ CD10 (56C6), 7 mL	PA0270	Bond™ Ready-to-Use Primary Antibody CD10 (56C6)
Bond™ CD11c (5D11), 7 mL	PA0554	Bond™ Ready-to-Use Primary Antibody CD11c (5D11)
Bond™ CD13 (38C12), 7 mL	PA0304	Bond™ Ready-to-Use Primary Antibody CD13 (38C12)
Bond™ CD15 (Carb-1), 7 mL	PA0039	Bond™ Ready-To-Use Primary Antibody CD15 (Carb-1)
Bond™ CD15 (MMA), 7 mL	PA0473	Bond™ Ready-to-Use Primary Antibody CD15 (MMA)
Bond™ CD19 (BT51E), 7 mL	PA0843	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		CD19 (BT51E)
Bond™ CD20 (L26), 7 mL	PA0200	Bond™ Ready-to-Use Primary Antibody CD20 (L26)
Bond™ CD20 (L26), 30 mL	PA0359	Bond™ Ready-to-Use Primary Antibody CD20 (L26)
Bond™ CD20 (MJ1), 7 mL	PA0906	Bond™ Ready-to-Use Primary Antibody CD20 (MJ1)
Bond™ CD21 (2G9), 7 mL	PA0171	Bond™ Ready-to-Use Primary Antibody CD21 (2G9)
Bond™ CD22 (FPC1), 7 mL	PA0249	Bond™ Ready-to-Use Primary Antibody CD22 (FPC1)
Bond™ CD23 (1B12), 7 mL	PA0169	Bond™ Ready-to-Use Primary Antibody CD23 (1B12)
Bond™ CD25 (4C9), 7 mL	PA0305	Bond™ Ready-to-Use Primary Antibody CD25 (4C9)
Bond™ CD30 (1G12), 7 mL	PA0153	Bond™ Ready-to-Use Primary Antibody CD30 (1G12)
Bond™ CD30 (JCM182), 7 mL	PA0790	Bond™ Ready-to-Use Primary Antibody CD30 (JCM182)
Bond™ CD31 (1A10), 7 mL	PA0250	Bond™ Ready-to-Use Primary Antibody CD31 (1A10)
Bond™ CD31 (JC70A), 7mL	PA0414	Bond™ Ready-to-Use Primary Antibody CD31 (JC70A)
Bond™ CD33 (PWS44), 7 mL	PA0555	Bond™ Ready-to-Use Primary Antibody CD33 (PWS44)
Bond™ CD34 (QBEnd/10), 7 mL	PA0212	Bond™ Ready-to-Use Primary Antibody CD34 (QBEnd/10)
Bond™ CD34 (QBEnd/10), 30 mL	PA0354	Bond™ Ready-to-Use Primary Antibody CD34 (QBEnd/10)
Bond™ CD43 (MT1), 7 mL	PA0938	Bond™ Ready-to-Use Primary Antibody CD43 (MT1)
Bond™ CD45 (X16/99), 7 mL	PA0042	Bond™ Ready-to-Use Primary Antibody CD45 (X16/99)
Bond™ CD45RO (UCHL1), 7 mL	PA0146	Bond™ Ready-To-Use Primary Antibody CD45RO (UCHL1)
Bond™ CD56 (CD564), 7 mL	PA0191	Bond™ Ready-to-Use Primary Antibody CD56 (CD564)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ CD57 (NK-1), 7 mL	PA0443	Bond™ Ready-To-Use Primary Antibody CD57 (NK-1)
Bond™ CD61 (2f2), 7 mL	PA0308	Bond™ Ready-to-Use Primary Antibody CD61 (2f2)
Bond™ CD68 (514H12), 7 mL	PA0273	Bond™ Ready-to-Use Primary Antibody CD68 (514H12)
Bond™ CD79a (11E3), 7 mL	PA0192	Bond™ Ready-to-Use Primary Antibody CD79a (11E3)
Bond™ CD79a (JCB117), 7 mL	PA0599	Bond™ Ready-to-Use Primary Antibody CD79a (JCB117)
Bond™ CD99 (12E7), 7 mL	PA0509	Bond™ Ready-To-Use Primary Antibody CD99 (12E7)
Bond™ CD103 (EP206), 7 mL	PA0374	Bond™ Ready-to-Use Primary Antibody CD103 (EP206)
Bond™ CD117 (EP10), 7mL	PA0007	Bond™ Ready-to-Use Primary Antibody CD117 (EP10)
Bond™ CD138 (Syndecan-1) (MI15), 7 mL	PA0088	Bond™ Ready-to-Use Primary Antibody CD138 (MI15)
Bond™ CD163 (10D6), 7 mL	PA0090	Bond™ Ready-to-Use Primary Antibody CD163 (10D6)
Bond™ CDX2 (AMT28), 7 mL	PA0535	Bond™ Ready-to-Use Primary Antibody CDX2 (AMT28)
Bond™ CDX2 (EP25), 7 mL	PA0375	Bond™ Ready-to-Use Primary Antibody CDX2 (EP25)
Bond™ Chromogranin A (5H7), 7 mL	PA0430	Bond™ Ready-to-Use Primary Antibody Chromogranin A (5H7)
Bond™ Chromogranin A (5H7), 7 mL	PA0515	Bond™ Ready-to-Use Primary Antibody Chromogranin A (5H7)
Bond™ Cyclin D1 (EP12), 7 mL	PA0046	Bond™ Ready-to-Use Primary Antibody Cyclin D1 (EP12)
Bond™ Cytokeratin 5 (XM26), 7 mL	PA0468	Bond™ Ready-to-Use Primary Antibody Cytokeratin 5 (XM26)
Bond™ Cytokeratin 7 (RN7), 7 mL	PA0942	Bond™ Ready-to-Use Primary Antibody Cytokeratin 7 (RN7)
Bond™ Cytokeratin 7 (RN7), 30mL	PA0138	Bond™ Ready-to-Use Primary Antibody Cytokeratin 7 (RN7)
Bond™ Cytokeratin 8 (TS1), 7 mL	PA0567	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Cytokeratin 8 (TS1)
Bond™ Cytokeratin 8/18 (5D3), 7 mL	PA0067	Bond™ Ready-to-Use Primary Antibody Cytokeratin 8/18 (5D3)
Bond™ Cytokeratin 14 (LL002), 7mL	PA0074	Bond™ Ready-to-Use Primary Antibody Cytokeratin 14 (LL002)
Bond™ Cytokeratin 17 (E3), 7mL	PA0114	Bond™ Ready-to-Use Primary Antibody Cytokeratin 17 (E3)
Bond™ Cytokeratin 19 (b170), 7 mL	PA0799	Bond™ Ready-to-Use Primary Antibody Cytokeratin 19 (b170)
Bond™ Cytokeratin 20 (Ks20.8), 7 mL	PA0022	Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (Ks20.8)
Bond™ Cytokeratin 20 (Ks20.8), 30 mL	PA0037	Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (Ks20.8)
Bond™ Cytokeratin 20 (PW31), 7 mL	PA0918	Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (PW31)
Bond™ Cytokeratin (AE1/AE3), 7mL	PA0094	Bond™ Ready-to-Use Primary Antibody Cytokeratin (AE1/AE3)
Bond™ Cytokeratin HMW (34βE12), 7 mL	PA0134	Bond™ Ready-to-Use Primary Antibody Cytokeratin HMW (34βE12)
Bond™ Desmin (DE-R-11), 7 mL	PA0032	Bond™ Ready-to-Use Primary Antibody Desmin (DE-R-11)
Bond™ DOG-1 (K9), 7 mL	PA0219	Bond™ Ready-to-Use Primary Antibody DOG-1 (K9)
Bond™ E-Cadherin (36B5), 7 mL	PA0387	Bond™ Ready-To-Use Primary Antibody E-Cadherin (36B5)
Bond™ Enhancer of Zeste Homolog 2 (Drosophila) (EZH2) (6A10), 7mL	PA0575	Bond™ Ready-to-Use Primary Antibody Enhancer of Zeste Homolog 2 (Drosophila) (EZH2) (6A10)
Bond™ Epithelial Membrane Antigen (GP1.4), 7 mL	PA0035	Bond™ Ready-to-Use Primary Antibody Epithelial Membrane Antigen (GP1.4)
Bond™ Estrogen Receptor (6F11), 7 mL	PA0151	Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody for Bond™
Bond™ Estrogen Receptor (6F11), 30 mL	PA0009	Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody for Bond™
Bond™ Factor XIIIa (E890.1), 7 mL	PA0449	Bond™ Ready-To-Use Primary Antibody Factor XIIIa (E980.1)
Bond™ Fascin (IM20), 7 mL	PA0420	Bond™ Ready-to-Use Primary Antibody Fascin (IM20)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ Galectin-3 (9C4), 7mL	PA0238	Bond™ Ready-to-Use Primary Antibody Galectin-3 (9C4)
Bond™ Gastrin (Polyclonal), 7mL	PA0681	Bond™ Ready-to-Use Primary Antibody Gastrin (Polyclonal)
Bond™ Glial Fibrillary Acidic Protein (GA5), 7 mL	PA0026	Bond™ Ready-to-Use Primary Antibody Glial Fibrillary Acidic Protein (GA5)
Bond™ Glucagon (Polyclonal), 7mL	PA0594	Bond™ Ready-to-Use Primary Antibody Glucagon (Polyclonal)
Bond™ Granzyme B (11F1), 7 mL	PA0291	Bond™ Ready-to-Use Primary Antibody Granzyme B (11F1)
Bond™ Gross Cystic Disease Fluid Protein-15 (23A3), 7 mL	PA0350	Bond™ Ready-to-Use Primary Antibody Gross Cystic Disease Fluid Protein-15 (23A3)
Bond™ Gross Cystic Disease Fluid Protein-15 (23A3), 7 mL	PA0708	Bond™ Ready-to-Use Primary Antibody Gross Cystic Disease Fluid Protein-15 (23A3)
Bond™ Human Chorionic Gonadotrophin (Polyclonal), 7 mL	PA0014	Bond™ Ready-To-Use Primary Antibody Human Chorionic Gonadotrophin (Polyclonal)
Bond™ Human Follicle Stimulating Hormone (INN-hFSH-60), 7 mL	PA0693	Bond™ Ready-to-Use Primary Antibody Human Follicle Stimulating Hormone (INN-hFSH-60)
Bond™ Human Growth Hormone (Polyclonal), 7 mL	PA0704	Bond™ Ready-to-Use Primary Antibody Human Growth Hormone (Polyclonal)
Bond™ Human Herpesvirus (Type 8) (13B10), 7mL	PA0050	Bond™ Ready-to-Use Primary Antibody Human Herpesvirus (Type 8)
Bond™ Immunoglobulin D (DRN1C), 7mL	PA0061	Bond™ Ready-to-Use Primary Antibody Immunoglobulin D (DRN1C)
Bond™ Immunoglobulin G (RWP49), 7mL	PA0905	Bond™ Ready-to-Use Primary Antibody Immunoglobulin G (RWP49)
Bond™ Immunoglobulin M (8H6), 7mL	PA0278	Bond™ Ready-to-Use Primary Antibody Immunoglobulin M (8H6)
Bond™ Inhibin Alpha (R1), 7 mL	PA0488	Bond™ Ready-to-Use Primary Antibody Inhibin Alpha (R1)
Bond™ Insulin (2D11-H5), 7 mL	PA0620	Bond™ Ready-to-Use Primary Antibody Insulin (2D11-H5)
Bond™ Ki67 (MM1), 7 mL	PA0118	Bond™ Ready-to-Use Primary Antibody Ki67 (MM1)
Bond™ Ki67 (MM1), 30 mL	PA0410	Bond™ Ready-to-Use Primary Antibody Ki67 (MM1)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ Ki67 (K2), 7 mL	PA0230	Bond™ Ready-to-Use Primary Antibody Ki67 (K2)
Bond™ Kappa Light Chain (CH15), 7mL	PA0606	Bond™ Ready-to-Use Primary Antibody Kappa Light Chain (CH15)
Bond™ Lambda Light Chain (SHL53), 7mL	PA0570	Bond™ Ready-to-Use Primary Antibody Lambda Light Chain (SHL53)
Bond™ Luteinizing Hormone (C93), 7 mL	PA0655	Bond™ Ready-to-Use Primary Antibody Luteinizing Hormone (C93)
Bond™ Macrophage Marker (MAC387), 7 mL	PA0752	Bond™ Ready-to-Use Primary Antibody Macrophage Marker (MAC387)
Bond™ Mammaglobin (EP249), 7 mL	PA0378	Bond™ Ready-to-Use Primary Antibody Mammaglobin (EP249)
Bond™ Mast Cell Tryptase (10D11), 7 mL	PA0019	Bond™ Ready-To-Use Primary Antibody Mast Cell Tryptase (10D11)
Bond™ Melan A (A103), 7 mL	PA0233	Bond™ Ready-to-Use Primary Antibody Melan A (A103)
Bond™ Melan A (A103), 30 mL	PA0044	Bond™ Ready-to-Use Primary Antibody Melan A (A103)
Bond™ Melanoma Marker (HMB45), 7mL	PA0027	Bond™ Ready-to-Use Primary Antibody Melanoma Marker (HMB45)
Bond™ Melanoma Marker (HMB45), 30 mL	PA0625	Bond™ Ready-to-Use Primary Antibody Melanoma Marker (HMB45)
Bond™ Mesothelin (5B2), 7 mL	PA0373	Bond™ Ready-to-Use Primary Antibody Mesothelin (5B2)
Bond™ MLH1 (Mismatch Repair Protein), (ES05), 7 mL	PA0610	Bond™ Ready-to-Use Primary Antibody MLH1 (Mismatch Repair Protein) (ES05)
Bond™ MSH2 (Mismatch Repair Protein) (25D12),7 mL	PA0048	Bond™ Ready-to-Use Primary Antibody MSH2 (Mismatch Repair Protein) (25D12)
Bond™ Muc-1 Glycoprotein (Ma695), 7mL	PA0051	Bond™ Ready-to-Use Primary Antibody Muc-1 Glycoprotein (Ma695)
Bond™ Muc-2 Glycoprotein (Ccp58), 7mL	PA0155	Bond™ Ready-to-Use Primary Antibody Muc-2 (Ccp58)
Bond™ Muc-5AC Glycoprotein (CLH2), 7mL	PA0052	Bond™ Ready-to-Use Primary Antibody Muc-5AC Glycoprotein (CLH2)
Bond™ Muc-6 Glycoprotein (CLH5), 7mL	PA0053	Bond™ Ready-to-Use Primary Antibody Muc-6 Glycoprotein (CLH5)
Bond™ Multi-Cytokeratin (AE1/AE3), 7 mL	PA0909	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Multi-Cytokeratin (AE1/AE3)
Bond™ Multi-Cytokeratin (AE1/AE3), 7 mL	PA0012	Bond™ Ready-to-Use Primary Antibody Multi-Cytokeratin (AE1/AE3)
Bond™ Multiple Myeloma Oncogene 1 (MUM1), (EAU32), 7 mL	PA0129	Bond™ Ready-to-Use Primary Antibody Multiple Myeloma Oncogene 1 (MUM1) (EAU32)
Bond™ Muramidase (Lysozyme) (Polyclonal), 7 mL	PA0391	Bond™ Ready-to-Use Primary Antibody Muramidase (Lysozyme) (Polyclonal)
Bond™ Muscle Specific Actin (HHF35), 7 mL	PA0258	Bond™ Ready-to-Use Primary Antibody Muscle Specific Actin (HHF35)
Bond™ Myeloperoxidase (59A5), 7 mL	PA0491	Bond™ Ready-to-Use Primary Antibody Myeloperoxidase (59A5)
Bond™ Myf-4 (Rhabdomyosarcoma Marker) (LO26), 7 mL	PA0226	Bond™ Ready-to-Use Primary Antibody Myf-4 (Rhabdomyosarcoma Marker) (LO26)
Bond™ Myoglobin (MYO18), 7 mL	PA0727	Bond™ Ready-to-Use Primary Antibody Myoglobin (MYO18)
Bond™ Myosin Heavy Chain (Smooth muscle) (S131), 7 mL	PA0493	Bond™ Ready-to-Use Primary Antibody Myosin Heavy Chain (smooth muscle) (S131)
Bond™ Napsin A (IP64), 7mL	PA0064	Bond™ Ready-to-Use Primary Antibody Napsin A (IP64)
Bond™ Negative (Rabbit), 7 mL	PA0777	Bond™ Ready-to-Use Negative Control Negative (Rabbit)
Bond™ Negative (Mouse), 7 mL	PA0996	Bond™ Ready-to-Use Negative Control Negative (Mouse)
Bond™ Neurofilament 200kD (N52.1.7), 7 mL	PA0371	Bond™ Ready-to-Use Primary Antibody Neurofilament 200 kD (N52.1.7)
Bond™ Neuron Specific Enolase (22C9), 7 mL	PA0435	Bond™ Ready-to-Use Primary Antibody Neuron Specific Enolase (22C9)
Bond™ Oct-2 (Oct-207), 7 mL	PA0532	Bond™ Ready-to-Use Primary Antibody Oct-2 (Oct-207)
Bond™ Oct-3/4 (N1NK), 7 mL	PA0934	Bond™ Ready-to-Use Primary Antibody Oct-3/4 (N1NK)
Bond™ Oct-3/4 (N1NK), 7 mL	PA0193	Bond™ Ready-to-Use Primary Antibody Oct-3/4 (N1NK)
Bond™ p120 Catenin (EP66), 7 mL	PA0379	Bond™ Ready-to-Use Primary Antibody p120 Catenin (EP66)

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
Bond™ p53 (DO-7), 7 mL	PA0057	Bond™ Ready-to-Use Primary Antibody p53 (DO-7)
Bond™ p63 (7JUL), 7 mL	PA0103	Bond™ Ready-to-Use Primary Antibody p63 (7JUL)
Bond™ Pax-5 (1EW), 7 mL	PA0552	Bond™ Ready-To-Use Primary Antibody Pax-5 (1EW)
Bond™ Placental Alkaline Phosphatase (8A9), 7 mL	PA0161	Bond™ Ready-to-Use Primary Antibody Placental Alkaline Phosphatase (8A9)
Bond™ Progesterone Receptor (16), 7 mL	PA0312	Bond™ Ready-to-Use Primary Antibody Progesterone Receptor (16)
Bond™ Programmed Death Ligand 1 (73-10), 7mL	PA0832	Bond™ Ready-to-Use Primary Antibody Programmed Death Ligand 1 (73-10)
Bond™ Prostatic Acid Phosphatase (PASE/4LJ), 7 mL	PA0006	Bond™ Ready-to-Use Primary Antibody Prostatic Acid Phosphatase (PASE/4LJ)
Bond™ Protein Gene Product 9.5 (10A1), 7 mL	PA0286	Bond™ Ready-to-Use Primary Antibody Protein Gene Product 9.5 (10A1)
Bond™ S-100 (Polyclonal), 7 mL	PA0900	Bond™ Ready-to-Use Primary Antibody S-100 (Polyclonal)
Bond™ Serotonin (Polyclonal), 7 mL	PA0736	Bond™ Ready-to-Use Primary Antibody Serotonin (Polyclonal)
Bond™ Smooth Muscle Actin (asm-1), 7 mL	PA0943	Bond™ Ready-to-Use Primary Antibody Smooth Muscle Actin (alpha sm-1)
Bond™ Somatostatin (Polyclonal), 7 mL	PA0331	Bond™ Ready-To-Use Primary Antibody Somatostatin (Polyclonal)
Bond™ Synaptophysin (27G12), 7 mL	PA0299	Bond™ Ready-to-Use Primary Antibody Synaptophysin (27G12)
Bond™ Tartrate-Resistance Acid Phosphatase (26E5), 7ml	PA0093	Bond™ Ready-to-Use Primary Antibody Tartrate-Resistant Acid Phosphatase (26E5)
Bond™ Terminal Deoxynucleotidyl Transferase (SEN28), 7 mL	PA0339	Bond™ Ready-to-Use Primary Antibody Terminal Deoxynucleotidyl Transferase (SEN28)
Bond™ Thyroglobulin (1D4), 7 mL	PA0025	Bond™ Ready-To-Use Primary Antibody Thyroglobulin (1D4)
Bond™ Thyroid Stimulating Hormone (QB2/6),7 mL	PA0776	Bond™ Ready-to-Use Primary Antibody Thyroid Stimulating Hormone (QB2/6)
Bond™ Thyroid Transcription Factor-1 (SPT24), 7 mL	PA0364	Bond™ Ready-to-Use Primary Antibody Thyroid Transcription Factor-1 (SPT24)
Bond™ Tyrosinase (T311), 7 mL	PA0322	Bond™ Ready-to-Use Primary Antibody

Device Schedule A (Continued)

Bond Ready To Use Antibodies	Product Code	IFU Product Name
		Tyrosinase (T311)
Bond™ Villin (CWWB1), 7 mL	PA0106	Bond™ Ready-to-Use Primary Antibody Villin (CWWB1)
Bond™ Vimentin (SRL33), 7 mL	PA0033	Bond™ Ready-To-Use Primary Antibody Vimentin (SRL33)
Bond™ Vimentin (V9) 7mL	PA0640	Bond™ Ready-to-Use Primary Antibody Vimentin (V9)
Bond™ von Willebrand Factor (36B11), 7 mL	PA0400	Bond™ Ready-to-Use Primary Antibody von Willebrand Factor (36B11)
Bond™ von Willebrand Factor (36B11), 7 mL	PA0055	Bond™ Ready-to-Use Primary Antibody von Willebrand Factor (36B11)
Bond™ Wilms' Tumor (WT49), 7 mL	PA0562	Bond™ Ready-to-Use Primary Antibody Wilms' Tumor (WT49)
Bond™ ZAP-70 (L453R), 7 mL	PA0998	Bond™ Ready-to-Use Primary Antibody ZAP-70 (L453R)

Consumables	Product Code
Bond™ Open Containers 7 ml	OP79193
Bond™ Open Containers 30 ml	OP309700
Bond™ Titration Kit	OPT9049
Bond™ Titration Container Inserts – 50 Pack	OPT9719

Signed: 

Date: 31 OCT 2018

Name: JEREMY ALLEN

Job Title: SNR QRA MANAGER

Letter of Authorization (LoA) for Commercialization and Regulatory Representation

To Whom It May Concern,

quartett Biotechnologie GmbH, with its registered address at Am Mühlenberg 4, 14476 Potsdam, Germany (hereinafter referred to as “**quartett**”), represented by its **CEO, Alexander Gorczyza**, hereby issues this **LoA**.

biocyc Biotechnologie GmbH & Co. KG, also located at Am Mühlenberg 4, 14476 Potsdam, Germany, and likewise represented by **CEO Alexander Gorczyza** (hereinafter “**biocyc**”), is acknowledged as the exclusive manufacturer of the products that quartett markets and distributes.

quartett hereby appoints **Medist Life Science S.R.L.**, located at **34, Ion Urdareanu Street, 050688, Sector 5, Bucharest, Romania**, to act on our behalf as a **non-exclusive commercial and regulatory representative** within the territory of **Romania and the Republic of Moldova**:


- Actively commercializing, promoting, and selling **quartett**'s medical device products within the specified territory,
- Assisting with registration and regulatory submissions related to our medical devices,
- Communicating with regulatory authorities and certification bodies,
- Supporting technical and administrative procedures,
- Representing us regarding compliance and conformity assessment

This authorization applies to all medical device products manufactured by **biocyc** and marketed by **quartett**, including—but not limited to—FFPE and frozen section antibodies, FITC-conjugated antibodies, immunohistochemistry reagents (buffers, antigen retrieval solutions, antibody diluents, mounting media, fixatives), as well as immunohistological detection kits.

This **LoA** becomes effective from **2026/01/02** and will remain valid until **2026/12/31**, unless earlier revoked by written notification from **quartett**.

This authorization does not transfer ownership rights or confer any exclusive rights regarding the marketing, distribution, or representation of our medical devices.

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declares under our sole responsibility that the products mentioned in this document have been manufactured and released according to the provisions of the QM system in force on the date of issue of this declaration and fulfil all relevant requirements of the following legal regulations and common specifications of the European Union:

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Anti-human antibody for *in vitro* diagnostic use. The primary antibody is intended for the qualitative detection of associated antigens as listed in the section ‘Summary and explanation’. It is intended to be used within an immunohistochemistry (IHC) procedure on formalin-fixed, paraffin-embedded (FFPE) tissue sections followed by light microscopy visualization to aid tumor diagnosis. The antibody may be used manually or with any automated staining platform.

Authorized and skilled personnel may only use the product. The clinical interpretation of any test results should be evaluated within the context of the patient’s medical history and other diagnostic laboratory test results. A qualified pathologist must perform evaluation.

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A handwritten signature in blue ink, appearing to read 'A. Gorczyza'.

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20ul selling size

Key facts

Isotype	IgG
Host species	Rabbit
Storage buffer	pH: 7.2 - 7.4 Preservative: 0.01% Sodium azide Constituents: PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA
Form	Liquid
Clonality	Monoclonal
Immunogen	The exact immunogen used to generate this antibody is proprietary information.
Clone number	EPR22811-18
Purification technique	Affinity purification Protein A
Concentration	0.51 mg/mL The concentration of this product may be batch-dependent Batch concentration finder →

Reactivity data

IP

Tested

Species	Human
Dilution info	1/30
Notes	-

Flow Cyt

Tested

Species	Human
Dilution info	1/500
Notes	-

WB

Tested

Species	Human
Dilution info	1/1000
Notes	-

ICC/IF

Tested

Species	Human
Dilution info	1/500
Notes	-

IHC-P

Tested

Species	Human
Dilution info	1/100
Notes	Perform heat-mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.

Target data

[See full target information ENG](#) 

Function Vascular endothelium glycoprotein that plays an important role in the regulation of angiogenesis (PubMed:21737454, PubMed:23300529). Required for normal structure and integrity of adult vasculature (PubMed:7894484). Regulates the

migration of vascular endothelial cells (PubMed:17540773). Required for normal extraembryonic angiogenesis and for embryonic heart development (By similarity). May regulate endothelial cell shape changes in response to blood flow, which drive vascular remodeling and establishment of normal vascular morphology during angiogenesis (By similarity). May play a critical role in the binding of endothelial cells to integrins and/or other RGD receptors (PubMed:1692830). Acts as a TGF-beta coreceptor and is involved in the TGF-beta/BMP signaling cascade that ultimately leads to the activation of SMAD transcription factors (PubMed:21737454, PubMed:22347366, PubMed:23300529, PubMed:8370410). Required for GDF2/BMP9 signaling through SMAD1 in endothelial cells and modulates TGFB1 signaling through SMAD3 (PubMed:21737454, PubMed:22347366, PubMed:23300529).

Storage

Shipped at conditions	Blue Ice
Appropriate short-term storage duration	1-2 weeks
Appropriate short-term storage conditions	+4°C
Appropriate long-term storage conditions	-20°C
Aliquoting information	Upon delivery aliquot
Storage information	Avoid freeze / thaw cycle

Notes

Patented technology

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.

What are the advantages of a recombinant monoclonal antibody?

This product is a recombinant monoclonal antibody, which offers several advantages including:

- High batch-to-batch consistency and reproducibility
- Improved sensitivity and specificity
- Long-term security of supply
- Animal-free batch production

For more information, read more on recombinant antibodies.

Supplementary info

This supplementary information is collated from multiple sources and compiled automatically.

Activity summary

CD105 also known as endoglin or the CD105 marker is a transmembrane glycoprotein with a molecular weight of approximately 180 kDa. It is a component of the TGF-beta receptor complex and exists in endothelial cells where it is abundantly expressed. Expression of CD105 is higher in proliferating cells particularly in the vasculature. You can also find it in tissues involved in the formation and remodeling of blood vessels such as during angiogenesis.

Biological function summary

Endoglin functions in the regulation of angiogenesis and vascular remodeling. It plays a significant role in mediating cellular responses to TGF-beta signaling influencing endothelial cell proliferation and migration. While not part of a larger structural complex endoglin interacts with receptors and signaling molecules important for vascular development and repair processes. This involvement aids in maintaining endothelial integrity and function under various physiological conditions.

Pathways

CD105 participates in the TGF-beta signaling and angiogenesis pathways. In these pathways it acts in conjunction with other proteins like TGF-beta receptors which play roles in cell differentiation proliferation and apoptosis. The interaction between CD105 and TGF-beta signaling regulates numerous cellular mechanisms impacting angiogenesis and cellular responses to environmental changes.

Associated diseases and disorders

CD105 has links to hereditary hemorrhagic telangiectasia (HHT) and certain cancers. In HHT mutations in the endoglin gene alter vascular structure leading to the formation of abnormal blood vessels. Oncologically overexpression of CD105 is present in tumor angiogenesis aiding in the progression of certain cancers. Other proteins like VEGF and TGF-beta closely interact with endoglin influencing disease progression and presenting potential targets for therapeutic intervention.

Product promise

Tested

We have tested this species and application combination and it works. It is covered by our product promise.

Expected

We have not tested this specific species and application combination in-house, but expect it will work. It is covered by our product promise.

Predicted

This species and application combination has not been tested, but we predict it will work based on strong homology. However, this combination is not covered by our product promise.

Not recommended

We do not recommend this combination. It is not covered by our product promise.


We are dedicated to supporting your work with high quality reagents and we are here for you every step of the way should you need us.

In the unlikely event of one of our products not working as expected, you are covered by our product promise.

Full details and terms and conditions can be found here:

Terms & Conditions.


8 product images

Flow Cytometry - Anti-CD105 antibody [EPR22811-18] (ab231774)

CD105 Flow Cytometry staining using rabbit Anti-CD105 antibody

Flow cytometric analysis of Jurkat (human T cell leukemia T lymphocyte, Left) / HUVEC (human umbilical vein endothelial cell, Right) cells labeling CD105 with ab231774 at 1/500 dilution (Red) compared with a Rabbit monoclonal IgG (ab172730) (Black) isotype control and an unlabeled control (cells without incubation with primary antibody and secondary antibody) (Blue). Goat anti rabbit IgG (Alexa Fluor® 488, ab150077) at 1/2000 dilution was used as the secondary antibody. **Negative control:** Jurkat (PMID: 28351936). Gated on viable cells.


Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD105 antibody [EPR22811-18] (ab231774)

CD105 Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) staining using rabbit Anti-CD105 antibody

Immunohistochemical analysis of paraffin-embedded human ovarian carcinoma tissue labeling CD105 with ab231774 at 1/100 dilution (5.1 µg/ml) followed by a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (ab209101). Positive staining on endothelial cells of human ovarian carcinoma (PMID: 17502949). The section was incubated with ab231774 for 30 mins at room temperature. The immunostaining was performed on a Leica Biosystems BOND® RX instrument. Counterstained with hematoxylin. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0, epitope retrieval solution 2) for 20mins.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (ab209101).


Immunocytochemistry/ Immunofluorescence - Anti-CD105 antibody [EPR22811-18] (ab231774)

CD105 Immunocytochemistry/ Immunofluorescence staining using rabbit Anti-CD105 antibody

Immunofluorescent analysis of 100% methanol-fixed U-937 (human histiocytic lymphoma monocyte) cells labeling CD105 with ab231774 at 1/500 dilution, followed by a ab150077 AlexaFluor®488 Goat anti-Rabbit secondary antibody at 1/1000 dilution (Green). Confocal image showing membranous staining in U-937 cells is observed. ab195889 Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) was used to counterstain tubulin at 1/200 dilution (Red). The nuclear counterstain was DAPI (Blue).

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is ab150077 AlexaFluor®488 Goat anti-Rabbit secondary at 1/1000 dilution.


Negative control: Jurkat (PMID: 28351936).

Flow Cytometry - Anti-CD105 antibody [EPR22811-18] (ab231774)

CD105 Flow Cytometry staining using rabbit Anti-CD105 antibody

Flow cytometric analysis of U-937 (human histiocytic lymphoma monocyte) cells labeling CD105 with ab231774 at 1/500 dilution (Red) compared with a Rabbit monoclonal IgG (ab172730) (Black) isotype control and an unlabeled control (cells without incubation with primary antibody and secondary antibody) (Blue). Goat anti rabbit IgG (Alexa Fluor® 488, ab150077) at 1/2000 dilution was used as the secondary antibody. Gated on viable cells.


Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD105 antibody [EPR22811-18] (ab231774)

CD105 Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) staining using rabbit Anti-CD105 antibody

Immunohistochemical analysis of paraffin-embedded human placenta tissue labeling CD105 with ab231774 at 1/100 dilution (5.1 µg/ml) followed by a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (ab209101). Positive staining on human placental trophoblasts (PMID: 17956952) is observed. The section was incubated with ab231774 for 30 mins at room temperature. The immunostaining was performed on a Leica Biosystems BOND® RX instrument. Counterstained with hematoxylin. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0, epitope retrieval solution 2) for 20mins.

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (ab209101).

Immunoprecipitation - Anti-CD105 antibody [EPR22811-18] (ab231774)

CD105 was immunoprecipitated from 0.35 mg HUVEC (human umbilical vein endothelial cell) whole cell lysate with ab231774 at 1/30 dilution (2µg in 0.35mg lysates). Western blot was performed on the immunoprecipitate using ab231774 1/1000 dilution (0.51 µg/ml). VeriBlot for IP secondary antibody (HRP) (ab131366) was used as the secondary antibody at 1/5000 dilution.

Lane 1: HUVEC (human umbilical vein endothelial cell) whole cell lysate 10µg

Lane 2: ab231774 IP in HUVEC whole cell lysate

Lane 3: Rabbit monoclonal IgG (ab172730) instead of ab231774 in HUVEC whole cell lysate.


Blocking and dilution buffer and concentration: 5% NFDN/TBST.

Exposure time: 3 mins.

All lanes:

Immunoprecipitation - Anti-CD105 antibody [EPR22811-18] (ab231774)

Predicted band size: 70 kDa

Western blot - Anti-CD105 antibody [EPR22811-18] (ab231774)

CD105 Western blot staining using rabbit Anti-CD105 antibody

Blocking/Diluting buffer and concentration: 5% NFDM/TBST.

Exposure times: Lanes 1-3: 26 secs; Lanes 4-5: 3 mins.

The molecular weight observed is consistent with what has been described in the literature (PMID: 8932339).

This blot was developed using a higher sensitivity ECL substrate.

Negative control: Jurkat (PMID: 28351936); Raji (PMID: 28351936).

All lanes:

Western blot - Anti-CD105 antibody [EPR22811-18] (ab231774) at 1/1000 dilution

Lane 1:

HUVEC (human umbilical vein endothelial cell), whole cell lysate 20 µg

Lane 2:

Jurkat (human T cell leukemia T lymphocyte), whole cell lysate 20 µg

Lane 3:

Raji (human Burkitt's lymphoma B lymphocyte), whole cell lysate 20 µg

Lane 4:

Human lung tissue lysate 20 µg

Lane 5:

Human placenta tissue lysate 20 µg

Secondary

Lanes 1 - 3:

Western blot - Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/100000 dilution

Lanes 4 - 5:

Western blot - VeriBlot for IP Detection Reagent (HRP) (ab131366) at 1/1000 dilution

Predicted band size: 70 kDa

Western blot - Anti-CD105 antibody [EPR22811-18] (ab231774)

CD105 Western blot staining of U937 (human histiocytic lymphoma monocyte), whole cell lysate 40 µg using rabbit Anti-CD105 antibody

Blocking/Diluting buffer and concentration: 5% NFDM/TBST.

The molecular weight observed is consistent with what has been described in the literature (PMID: 8932339).

All lanes:

Western blot - Anti-CD105 antibody [EPR22811-18] (ab231774) at 1/1000 dilution

All lanes:

U937 (human histiocytic lymphoma monocyte), whole cell lysate 40 µg

Secondary

All lanes:

Western blot - Goat Anti-Rabbit IgG H&L (HRP) (ab97051) at 1/20000 dilution

Predicted band size: 70 kDa

Observed band size: 97 kDa

Exposure time: 103s

Please note: All products are 'FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES'.

PAX-8 (MRQ-50) Mouse Monoclonal Antibody

Product Availability

REF	Description
PA0808	PAX8 (MRQ-50 Mab) 'CM', 7 mL

Symbol Definitions

KEY-CODE	keycode
Rab	rabbit antibody
Mab	mouse antibody

Intended Use

PAX-8 (MRQ-50) Mouse Monoclonal Primary Antibody is intended for laboratory use in the detection of the PAX-8 protein in formalin-fixed, paraffin-embedded tissue stained on Leica Bond IHC instruments. This product should be interpreted by a qualified pathologist in conjunction with histological examination, relevant clinical information, and proper controls. This antibody is intended for *in vitro* diagnostic (IVD) use.

Summary and Explanation

PAX-8 is a transcription factor expressed during embryonic development of Müllerian organs, kidney, and thyroid, with continued expression in some epithelial cell types of these mature tissues.¹ It can be useful for marking several types of carcinoma including ovarian serous carcinoma, clear cell renal cell carcinoma, and papillary thyroid carcinoma.¹⁻⁵ Additionally, PAX-8 is not found in the epithelial cells of the breast, lung, mesothelium, stomach, colon, pancreas and other sites.¹⁻⁴

Principles and Procedures

The stated primary antibody (PA0808) may be used as the primary antibody for immunohistochemical staining of formalin-fixed, paraffin-embedded tissue sections. In general, immunohistochemical staining in conjunction with a streptavidin-biotin detection system allows the visualization of antigens via the sequential application of a specific antibody (primary antibody) to the antigen, a secondary antibody (link antibody) to the primary antibody, an enzyme complex and a chromogenic substrate with interposed washing steps. Alternatively, a biotin-free detection system may be used. The enzymatic activation of the chromogen results in a visible reaction product at the antigen site. The specimen may then be counterstained and a coverslip applied. Results are interpreted using a light microscope and aid in the differential diagnosis of pathophysiological processes, which may or may not be associated with a particular antigen.

The stated primary antibody (PA0808) is optimally diluted to be compatible with Leica BOND System detection kits and automated slide stainers. For more detailed information on instrument operation, refer to the appropriate Leica BOND automated slide stainer Operator's Manual.

Materials and Methods

Reagents Provided

One dispenser of the stated primary antibody (PA0808) contains sufficient prediluted reagent for 45 tests. This antibody is diluted in Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide.

The immunoglobulin concentration range of the reagent appears on the product label.

Product Composition	
Predilute: diluted in	Tris Buffer, pH 7.3-7.7, with 1% BSA and <0.1% Sodium Azide
Host	Mouse
Isotype	IgG
Source	Supernatant

Reconstitution, Mixing, Dilution, Titration

This antibody is optimized for use on a Leica BOND automated slide stainer in combination with Leica BOND detection systems. No reconstitution, mixing, dilution, or titration is required. Further dilution may result in loss of antigen staining. The user must validate any such changes. Differences in tissue processing and technical procedures in the laboratory may produce significant variability in results and require regular use of controls. (See Quality Control Procedures section)

Materials and Reagents Needed But Not Provided

The following reagents and materials may be required for staining but are not provided with the primary antibody:

1. Positive and negative control tissue
2. Microscope slides, positively charged
3. Drying oven capable of maintaining a temperature of 58-60°C ± 5° C
4. Bar code labels (appropriate bar code labels for negative control and the primary antibody being tested)
5. Staining jars or baths
6. Timer
7. Amplifier (when applicable)
8. Xylene or xylene substitute
9. Ethanol or reagent alcohol
10. Deionized or distilled water
11. Leica BOND-MAX and Leica BOND-III autostainers
12. See the "Using BOND Reagents" section in the Leica BOND user documentation for complete listing of materials for use with BOND autostainers.
13. Hematoxylin or other counterstain
14. Negative Control Reagent

15. Mounting medium
16. Cover glass
17. Light microscope (40–400x)

Storage and Handling

Store at 2–8° C. Do not freeze.

To ensure proper reagent delivery and stability of the antibody after every run, the cap must be replaced and the dispenser must be immediately placed in the refrigerator in an upright position.

Every antibody dispenser is expiration dated. When properly stored, the reagent is stable to the date indicated on the label. Do not use reagent beyond the expiration date for the prescribed storage method.

There are no definitive signs to indicate instability of this product; therefore, positive and negative controls should be run simultaneously with unknown specimens. Contact Cell Marque technical support if there is a suspected indication of reagent instability.

Specimen Collection and Preparation for Analysis

Routinely processed, neutral-buffered formalin-fixed, paraffin-embedded tissues are suitable for use with this primary antibody when used with BOND Polymer Refine detection systems and BOND ancillary reagents on the Leica BOND Platforms (see Materials and Reagents Needed, But Not Provided section). The recommended tissue fixative is 10% neutral buffered formalin. Variable results may occur as a result of prolonged fixation or special processes such as decalcification of bone marrow preparations.

Each section should be cut to the appropriate thickness (approximately 3 µm) and placed on a positively charged glass slide. Slides containing the tissue section may be baked for at least 2 hours (but not longer than 24 hours) in a 58–60° C ± 5° C oven.

Warnings and Precautions

1. Take reasonable precautions when handling reagents. Use disposable gloves and lab coats when handling suspected carcinogens or toxic materials (example: xylene).
2. Avoid contact of reagents with eyes and mucous membranes. If reagents come in contact with sensitive areas, wash with copious amounts of water.
3. Patient specimens and all materials contacting them should be handled as biohazardous materials and disposed of with proper precautions. Never pipette by mouth.
4. Avoid microbial contamination of reagents, as this could produce incorrect results.
5. Incubation times and temperatures other than those specified may give erroneous results.
6. The reagents have been optimally diluted, and further dilution may result in loss of antigen staining. The user must validate any such change.
7. When used according to instructions, this product is not classified as a hazardous substance. The preservative in the reagent is less than 0.1% sodium azide and does not meet the OSHA (USA) criteria for hazardous substance at the stated concentration. See MSDS.
8. The user must validate any storage conditions other than those specified in the package insert.
9. Diluent may contain bovine serum albumin and supernatant may contain bovine serum. The products containing fetal bovine serum and products containing bovine serum albumin are purchased from commercial suppliers. Certificates of Origin for the animal source used in these products are on file at Cell Marque. The certificates support that the

bovine sources are from countries with negligible BSE risk and state sources of bovine from USA and Canada.

10. As with any product derived from biological sources, proper handling procedures should be used.

Instructions For Use

Dispenser Storage & Handling

To insure reliable operation, the dispenser must always be capped when not in use and should NEVER be manually dispensed.

Recommended Staining Protocols for PA0808

PAX-8 (MRQ-50) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection.

HRP Detection:

The recommended HRP staining protocol for PAX-8 (MRQ-50) primary antibody is IHC Protocol F.

Heat induced epitope retrieval is recommended using BOND™ Epitope Retrieval Solution 2 for 10 minutes.

Quality Control Procedures

Positive Tissue Control

A positive tissue control must be run with every staining procedure performed. This tissue may contain both positive and negative staining cells or tissue components and serve as both the positive and negative control tissue. Control tissues should be fresh autopsy, biopsy or surgical specimens prepared or fixed as soon as possible in a manner identical to the test sections. Use of a tissue section fixed or processed differently from the test specimen will serve to provide control for all reagents and method steps except fixation and tissue processing.

A tissue with weak positive staining is more suitable for optimal quality control and for detecting minor levels of reagent degradation. Positive tissue control for the stated primary antibody (PA0808) may include the following:

Positive Tissue Control	
Tissue	Visualization
Ovarian carcinoma (non-mucinous carcinoma)	Nuclear
Thyroid carcinoma	Nuclear
Renal cell carcinoma	Nuclear

Known positive tissue controls should be utilized only for monitoring the correct performance of processed tissues and test reagents, not as an aid in determining a specific diagnosis of patient samples. If the positive tissue controls fail to demonstrate appropriate positive staining, results with the test specimens must be considered invalid.

Negative Tissue Control

The same tissue used for the positive tissue control may be used as the negative tissue control. The variety of cell types present in most tissue sections offers internal negative control sites, but this should be verified by the user. The components that do not stain should demonstrate the absence of specific staining, and provide an indication of non-specific background staining. If specific staining occurs in the negative tissue control sites, results with the patient specimens must be considered invalid.

Unexplained Discrepancies

Unexplained discrepancies in controls should be referred to Leica Biosystems. If quality control results do not meet specifications, patient results are invalid. See the Troubleshooting section of this insert. Identify and correct the problem, then repeat the entire procedure with the patient samples.

Negative Control Reagent

A negative control reagent must be run for every specimen to aid in the interpretation of results. A negative control reagent is used in place of the primary antibody to evaluate nonspecific staining. The slide should be treated with negative control reagent, matching the host species of the primary antibody, and ideally having the same IgG concentration. The incubation period for the negative control reagent should equal the primary antibody incubation period.

Interpretation of Results

The immunostaining procedure run on Leica BOND automated slide stainers causes a colored reaction product to precipitate at the antigen sites localized by the primary antibody. Refer to the appropriate detection system package insert for expected color reactions. A qualified pathologist experienced in immunohistochemistry procedures must evaluate positive and negative tissue controls before interpreting results.

Positive Tissue Control

The stained positive tissue control should be examined first to ascertain that all reagents are functioning properly. The presence of an appropriately colored reaction product within the target cells is indicative of positive reactivity. Refer to the package insert of the detection system used for expected color reactions. Depending on the incubation length and potency of the hematoxylin used, counterstaining will result in a pale to dark blue coloration of cell nuclei. Excessive or incomplete counterstaining may compromise proper interpretation of results. If the positive tissue control fails to demonstrate appropriate positive staining, any results with the test specimens are considered invalid.

Negative Tissue Control

The negative tissue control should be examined after the positive tissue control to verify the specific labeling of the target antigen by the primary antibody. The absence of specific staining in the negative tissue control confirms the lack of antibody cross reactivity to cells or cellular components. If specific staining occurs in the negative tissue control, results with the patient specimen are considered invalid. Nonspecific staining, if present, will have a diffuse appearance. Sporadic light staining of connective tissue may also be observed in sections from tissues that are not optimally fixed. Intact cells should be used for interpretation of staining results. Necrotic or degenerated cells show non-specific staining.

Patient Tissue

Patient specimens should be examined last. Positive staining intensity should be assessed within the context of any background staining of the negative reagent control. As with any

immunohistochemical test, a negative result means that the antigen in question was not detected, not that the antigen is absent in the cells or tissue assayed. A panel of antibodies may aid in the identification of false negative reactions (see Summary of Expected Results section). The morphology of each tissue sample should also be examined utilizing a hematoxylin and eosin stained section when interpreting any immunohistochemical result. The patient's morphologic findings and pertinent clinical data must be interpreted by a qualified pathologist.

Limitations

1. This reagent is "for professional use only" as immunohistochemistry is a multiple step process that requires specialized training in the selection of the appropriate reagents, tissues, fixation, processing; preparation of the immunohistochemistry slide; and interpretation of the staining results.
2. For laboratory use only.
3. For *in vitro* diagnostic use.
4. Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning, or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may result from variations in fixation and embedding methods, as well as from inherent irregularities within the tissue.
5. Excessive or incomplete counterstaining may compromise proper interpretation of results.
6. The clinical interpretation of any positive staining, or its absence, must be evaluated within the context of clinical history, morphology, other histopathological criteria as well as other diagnostic tests. This antibody is intended to be used in a panel of antibodies. It is the responsibility of a qualified pathologist to be familiar with the antibodies, reagents and methods used to produce the stained preparation. Staining must be performed in a certified, licensed laboratory under the supervision of a pathologist who is responsible for reviewing the stained slides and assuring the adequacy of positive and negative controls.
7. Cell Marque provides antibodies at optimal dilution for use as instructed. Any deviation from recommended test procedures may invalidate expected results. Appropriate controls must be employed and documented. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results.
8. This product is not intended for use in flow cytometry; performance characteristics have not been determined.
9. Reagents may demonstrate unexpected reactions in previously untested tissues. The possibility of unexpected reactions even in tested tissue groups cannot be completely eliminated because of biological variability of antigen expression in neoplasms, or other pathological tissues. Contact Cell Marque technical support with documented unexpected reactions.
10. Tissues from persons infected with hepatitis B virus and containing hepatitis B surface antigen (HBsAg) may exhibit nonspecific staining with horseradish peroxidase.
11. When used in blocking steps, normal sera from the same animal source as the secondary antisera may cause false negative or false positive results because of the effect of autoantibodies or natural antibodies.
12. False positive results may be seen because of nonimmunological binding of proteins or substrate reaction products. They may also be caused by pseudoperoxidase activity (erythrocytes), endogenous peroxidase activity (cytochrome C), or endogenous biotin (example: liver, brain, breast, kidney) subject to the type of immunostaining technique used.

13. As with any immunohistochemistry test, a negative result means that the antigen was not detected, not that the antigen was absent in the cells or tissue assayed.

Specific Limitations

- The antibody is optimized for the incubation time specified in the Instructions for Use section in combination with Leica BOND System detection kits and the Leica BOND automated slide stainers. Because of variation in tissue fixation and processing, it may be necessary to increase or decrease the primary antibody incubation time on individual specimens.
- Cell Marque antibodies, when used in combination with Leica BOND detection systems and accessories, detects antigen(s) that survive routine formalin fixation, tissue processing, and sectioning. Users who deviate from recommended test procedures are responsible for interpretation and validation of patient results.

Summary of Expected Results

See the following tables of reactivity:

Normal Study			
Tissue	# Positive Cases	Total Cases Tested	Notes
Brain	0	1	
Adrenal Cortex	0	1	
Ovary	0	1	
Pancreas	1	1	Islet cells +
Parathyroid	0	1	
Pituitary	0	1	
Testis	0	1	
Thyroid	1	1	
Breast	0	1	
Spleen	1	1	
Tonsil	1	1	
Thymus	1	1	Mature lymphocytes +
Bone Marrow	0	1	
Lung	0	1	
Heart	0	1	
Esophagus	0	1	
Stomach	0	1	
Small Intestine	0	1	
Colon	0	1	
Liver	0	1	
Salivary Gland	0	1	
Gall Bladder	0	1	

Normal Study			
Tissue	# Positive Cases	Total Cases Tested	Notes
Kidney	1	1	
Bladder	0	1	
Prostate	0	1	
Uterus	1	1	
Fallopian Tube	1	1	
Ureter	1	1	Squamous epithelium +
Cervix	0	1	
Skeletal Muscle	0	1	
Smooth Muscle	0	1	
Skin	0	1	
Peripheral Nerve	0	1	
Mesothelium	0	1	
Fat	0	1	
Placenta	0	1	

This antibody stains normal tissues as indicated in literature.

Disease Tissue Study			
Tissue	# Positive Cases	Total Cases Tested	Notes
Renal cell carcinoma	26	27	
Ovarian carcinoma	18	20	
Papillary thyroid carcinoma	5	9	
Poorly differentiated thyroid carcinoma	1	1	
Breast carcinoma	0	5	
Breast invasive ductal carcinoma	0	22	
Chromophobe renal cell carcinoma	0	2	
Colorectal carcinoma	0	20	
Lung adenocarcinoma	0	5	
Transitional cell carcinoma	0	1	

This antibody stains tumors as indicated in literature.

Troubleshooting

1. If the positive control exhibits weaker staining than expected, other positive controls run during the same instrument run should be checked to determine if it is because of the primary antibody or one of the common secondary reagents.
2. If the positive control is negative, it should be checked to ensure that the slide has the proper bar code label. If the slide is labeled properly, other positive controls used on the same instrument run should be checked to determine if it is because of the primary antibody or one of the common secondary reagents. Tissues may have been improperly collected, fixed or deparaffinized. The proper procedure should be followed for collection, storage and fixation.
3. If excessive background staining occurs, high levels of endogenous biotin may be present. A biotin blocking step should be included unless a biotin-free detection system is being used in which case any biotin present would not be a contributing factor to background staining.
4. If all of the paraffin has not been removed, the deparaffinization procedure should be repeated.
5. If specific antibody staining is too intense, the run should be repeated with incubation time shortened by 4 minute intervals until the desired stain intensity is achieved.
6. If tissue sections wash off the slide, slides should be checked to ensure that they are positively charged.

For corrective action, refer to the Step By Step Procedure section, the automated slide stainer Operator's Manual, or contact Leica Biosystems customer service.

References

- 1 Ozcan A, et al. PAX 8 expression in non-neoplastic tissues, primary tumors, and metastatic tumors: a comprehensive immunohistochemical study. *Mod Pathol*. 2011; 24:751-64.
- 2 Laury AR, et al. A comprehensive analysis of PAX8 expression in human epithelial tumors. *Am J Surg Pathol*. 2011; 35:816-26.
- 3 Nonaka, D et al. Expression of pax8 as a useful marker in distinguishing ovarian carcinomas from mammary carcinomas. *Am J Surg Pathol*. 2008; 32:1566-71.
- 4 Nonaka, D et al. Diagnostic utility of thyroid transcription factors Pax8 and TTF-2 (Fox E1) in thyroid epithelial neoplasms. *Mod Pathol* .2008; 21:192-200.
- 5 Tong, GX et al. Expression of PAX8 in normal and neoplastic renal tissues: an immunohistochemical study. *Mod Pathol*. 2009; 22:1218-27.

Disclaimers

* Leica BOND-MAX system and Leica BOND-III are trademarks of Leica Biosystems. Cell Marque antibodies are developed manufactured and distributed by Cell Marque Corporation and their sale through Leica Microsystems does not imply approval, endorsement, or any guarantee of quality or performance of those Cell Marque antibodies by Leica Biosystems.

 www.cellmarque.com

 6600 Sierra College Blvd. • Rocklin, CA 95677 USA • 916-746-8900

CM Template #0.1
Implementation date 22 Jan 2018

BOND Ready-to-Use Primary Antibody Wilms' Tumor (WT49)

Catalog No: PA0562

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
+44 191 215 4242

[EN](#) [FR](#) [IT](#) [DE](#) [ES](#) [PT](#) [SV](#) [EL](#) [DA](#) [NL](#)

[NO](#) [TR](#) [BG](#) [HU](#) [RO](#) [RU](#) [PL](#) [SL](#) [CS](#) [SK](#)

[AR](#) [SR](#) [LV](#) [LT](#) [ET](#)

Instructions for Use

Please read before using this product.

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Catalog No: PA0562

Intended Use

This reagent is for *in vitro* diagnostic use.

Wilms' Tumor (WT49) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human Wilms' Tumor WT1 gene product in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes BOND-MAX system, BOND-III system and BOND-PRIME system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Wilms' Tumor (WT49) primary antibody is a ready to use product that has been optimized for use with the BOND Polymer Refine Detection system (DS9800) and BOND-PRIME Polymer DAB Detection System (DS9284). The demonstration of human Wilms' Tumor WT1 gene product is achieved by first allowing the binding of Wilms' Tumor (WT49) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes BOND-MAX system, BOND-III system and BOND-PRIME system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Wilms' Tumor (WT49) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

WT49.

Immunogen

A prokaryotic recombinant protein containing 1-181 amino acids of the N-terminal of the Wilms' Tumor protein.

Specificity

Human Wilms' Tumor WT1 gene product.

Subclass

IgG1.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 2.3 mg/L as determined by ELISA.

Dilution and Mixing

Wilms' Tumor (WT49) primary antibody is optimally diluted for use on the BOND system (includes BOND-MAX system, BOND-III system and BOND-PRIME system). Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation, or Sections 1 and 3 in your BOND-PRIME user documentation, for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes BOND-MAX system, BOND-III system and BOND-PRIME system).

Storage and Stability

Store at 2-8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Wilms' Tumor (WT49) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2-8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, LeicaBiosystems.com.

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Wilms' Tumor (WT49) primary antibody was developed for use on the automated BOND system (includes BOND-MAX system, BOND-III system and BOND-PRIME system) in combination with BOND Polymer Refine Detection and BOND-PRIME Polymer DAB Detection System. The recommended staining protocol and epitope retrieval for Wilms' Tumor (WT49) primary antibody are detailed in Table 1.

Table 1: Protocol Parameters for each BOND System.

	BOND-MAX	BOND-III	BOND-PRIME
Epitope Retrieval	Heat induced epitope retrieval using BOND Epitope Retrieval Solution 2 for 30 minutes	Heat induced epitope retrieval using BOND Epitope Retrieval Solution 2 for 30 minutes	Heat induced epitope retrieval using BOND-PRIME Epitope Retrieval Solution 2 for 30 minutes
Staining Protocol	*IHC Protocol F	*IHC Protocol F	*IHC Protocol F

Results Expected (generated with *IHC Protocol F on BOND-III platform using BOND Polymer Refine Detection System)

Normal Tissues

Clone WT49 detected the Wilms' tumor antigen in the nuclei of stromal cells in ovary, uterus, and cervix, in Sertoli cells of the testis and in podocytes and glomeruli of kidney. Clone WT49 does not stain endothelial cells. (Total number of cases evaluated = 139).

Tumor Tissues

Clone WT49 stained 4/85 kidney tumors (including 4/9 Wilms' tumors, 0/32 clear cell carcinomas, 0/18 transitional cell carcinomas, 0/10 papillary renal cell carcinomas, 0/5 squamous cell carcinomas, 0/2 collecting duct carcinomas, 0/2 undifferentiated carcinomas, 0/4 metastatic clear cell carcinomas, 0/1 metastatic chromophobe carcinoma, and 0/2 hyperplasias), 4/5 renal inflammatory conditions (including 2/3 interstitial nephritis, 1/1 chronic pyelonephritis, and 1/1 acute pyelonephritis), and 5/7 ovarian tumors (including 4/4 serous tumors, 1/1 endometrioid adenocarcinoma, 0/1 adenocarcinoma and 0/1 granulosa cell tumor). No staining was seen in lung tumors (0/5), thyroid tumors (0/5), carcinomas of the cervix tumors (0/4), colon tumors (0/4), hepatocellular carcinomas (0/4), breast tumors (0/4), brain tumors (0/4), lymphomas (0/3), adrenal tumors (0/2), bladder tumors (0/2), esophageal tumors (0/2), gastric tumors (0/2), tumors of the small intestine (0/2), rectal tumors (0/3), metastatic tumors (0/3), head and neck tumors (0/2), prostatic tumors (0/2), salivary gland tumors (0/2), seminomas (0/2), a chondrosarcoma (0/1), a melanoma (0/1), a tongue tumor (0/1), a pancreatic tumor (0/1), a skin tumor (0/1), and a hyperplastic prostate (0/1). (Total number of cases evaluated = 160).

Wilms' Tumor (WT49) is recommended for the detection of the Wilms' tumor protein in normal and tumor tissues

Product Specific Limitations

Wilms' Tumor (WT49) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection, BOND ancillary reagents, BOND-PRIME Polymer DAB Detection System, and BOND-PRIME ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProCin™ 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

Change History

Revision: Date of Issue	Detail of Revision
17 October 2022 (A)	<p>Change History: Inclusion of new section containing Change History table and Revision: Date of Issue.</p> <p>Revision / Date of Issue: Change of title of section from Date of Issue.</p> <p>General Changes: Removal of www. from LeicaBiosystems.com.</p> <p>Results Expected: Updated with evaluated tissues and data.</p> <p>Front Page: Removal of Rx Only.</p> <p>Back Page: Addition of EC REP details and symbol.</p>
B	US address updated. US and Canada telephone number updated. Australian telephone number removed. Name and address of EU Rep updated, correction in name of Canada entity.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 est une marque commerciale de Supelco, qui est une filiale de Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 è un marchio registrato di Supelco, parte di Sigma-Aldrich Corporation.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 ist eine eingetragene Marke von Supelco, einem Unternehmen der Sigma-Aldrich Corporation.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 es una marca comercial de Supelco, parte de Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 är ett varumärke som tillhör Supelco, en del av Sigma-Aldrich Corporation.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
Το ProClip™ 950 αποτελεί εμπορικό σήμα της Supelco, ενός τμήματος της Sigma-Aldrich Corporation.

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Litteraturliste

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 er et varemerke for Supelco, en del av Sigma-Aldrich Corporation.

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Kaynakça

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950, Sigma-Aldrich Corporation'a bağlı Supelco'nun ticari markasıdır.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 е търговска марка на Supelco, част от Sigma-Aldrich Corporation.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- A ProCin™ 950 a Sigma-Aldrich Corporation részét képező Supelco védjegye.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 — торговая марка компании Supelco, входящей в состав корпорации Sigma-Aldrich.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 jest znakiem handlowym firmy Supelco, będącej częścią Sigma-Aldrich Corporation.

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Literatura

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProCin™ 950 je blagovna znamka družbe Supelco, ki je del družbe Sigma-Aldrich Corporation.

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Literatúra

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 je ochranná známka spoločnosti Supelco, súčasti spoločnosti Sigma-Aldrich Corporation.

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Bibliografija

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Omeroglu A and Omeroglu G. Pathologic quiz case: a 13-year-old with multiple abdominal masses. Arch Pathol Lab Med. 2003; 127:e347-e348
 5. Lee SB, Haber DA. Wilms tumor and the WT1 gene. Experimental Cell Research. 2001;264(1):74-99.
- ProClin™ 950 jeste zaštitni znak kompanije Supelco, koja je deo korporacije Sigma-Aldrich.

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Leica Microsystems Canada
71 Four Valley Drive
Concord, Ontario L4K 4V8
Canada
J +1 844 534 2262

Leica Biosystems, division of Leica Microsystems, Inc.
21440 W. Lake Cook Road
Floor 5
Deer Park, IL 60010
USA
J +1 844 534 2262

Leica Biosystems Melbourne Pty Ltd
495 Blackburn Road
Mt Waverley VIC 3149
Australia

AR Experts B.V.
Boeingavenue 209
1119 PD, Schiphol-Rijk
The Netherlands

PAX8 (QR016)

EN	INSTRUCTIONS FOR USE	2
DE	GEBRAUCHSANWEISUNG	6
FR	MODE D'EMPLOI	10
ES	INSTRUCCIONES DE USO	14
IT	ISTRUZIONI PER L'USO	18
PT	INSTRUÇÕES DE UTILIZAÇÃO	22
GR	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	26
NO	BRUKSANVISNING	30
SV	BRUKSANVISNING	34
DA	BRUGSANVISNING	38
ET	KASUTUSJUHEND	42
LV	LIETOŠANAS INSTRUKCIJA	46
LT	NAUDOJIMO INSTRUKCIJOS	50
CS	NÁVOD K POUŽITÍ	54
HU	HASZNÁLATI UTASÍTÁS	58
RO	INSTRUCȚIUNI DE UTILIZARE	62
SK	NÁVOD NA POUŽITIE	66
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА	70
HR	UPUTE ZA UPORABU	74
TR	KULLANIM TALIMATLARI	78

EN INSTRUCTIONS FOR USE

Rabbit monoclonal antibody against PAX8 (QR016)

In Vitro Diagnostic Use (IVD)


Figure 1 Ovarian carcinoma stained with anti-PAX8 (QR016)

Product identification

C-P008-025	25 µl	Concentrate
C-P008-01	0.1 ml	Concentrate
C-P008-05	0.5 ml	Concentrate
C-P008-10	1 ml	Concentrate
P-P008-30	3 ml	Ready-to-use
P-P008-70	7 ml	Ready-to-use
P-P008-150	15 ml	Ready-to-use

Intended use

Anti-human antibody for *in vitro* diagnostic use. The primary antibody is intended for the qualitative detection of associated antigens as listed in the section 'Summary and explanation'. It is intended to be used within an immunohistochemistry (IHC) procedure on formalin-fixed, paraffin-embedded (FFPE) tissue sections followed by light microscopy visualization to aid tumor diagnosis. The antibody may be used manually or with any automated staining platform.

Authorized and skilled personnel may only use the product. The clinical interpretation of any test results should be evaluated within the context of the patient's medical history and other diagnostic laboratory test results. A qualified pathologist must perform evaluation.

Summary and explanation

PAX8 is a member of the PAX family (Paired Box) of transcription factors. It is vital for organogenesis during the embryonic development of the kidneys, Müller organs and thyroid. Due to the restrictive expression in normal tissues, PAX8 is a sensitive and specific marker for primary tumors as well as for metastatic tumors from the above-mentioned organs and tissues.

PAX8 is expressed in thyroid carcinoma (~90%), endometrial carcinoma (84-98%), ovarian carcinoma (71-99%) and renal cell carcinoma (~90%)^[1-2].

Principle of the procedure

The stated primary antibody is suitable for immunohistochemical staining of FFPE tissue sections based on specific antigen-antibody reaction. Using a detection system linked to horseradish peroxidase (HRP) or alkaline phosphatase (AP) the antigen visualization is performed via specific binding of the primary antibody.

Secondary antibody is binding to the primary antibody, and the enzyme complex labels this complex. The enzymatic activation of the chromogen results in a visible reaction product at the antigen site. Each step is incubated for a precise time and temperature and requires interposed washing steps. The specimen may then be counterstained. Results are interpreted using a light microscope.

Materials provided

Primary antibody	Anti-PAX8 (QR016)
Host	Rabbit
Subclass	IgG
Immunogen	Synthetic peptide of human PAX8
Antibody concentrate	Concentrated antibody in TRIS (pH 7.4) with < 1 % sera (bovine, donkey) and < 0.1 % sodium azide
Recommended working dilution range	1:100 – 1:200
Ready-to-use antibody	Prediluted antibody in TRIS (pH 7.4) with < 1 % sera (bovine, donkey) and < 0.1 % sodium azide

Product label shows the specific lot number. Each individual lot is compared and adjusted to a reference lot to ensure a consistent immunohistochemical staining performance from lot-to-lot.

Prediluted antibody is ready-to-use and optimized for staining. No further dilution, reconstitution, mixing, or titration is needed.

Antibody concentrate is optimized for dilution within dilution range using Q Diluent for IHC (Cat. No. AD-001-xxxx). Indicated dilution range should be considered as recommendation and depends on different facts (tissue, fixation, incubation conditions, etc.). Optimum dilution to be determined in user's own system.

Materials required but not provided

- Positive and negative controls
- Microscope slides (positively charged) and cover slips
- Staining jars
- Timer
- Xylene or xylene alternative, e.g. Q Dewax Solution (Cat. No. DW-001-xxxx)
- Ethanol
- Deionized or distilled water
- Heating equipment for tissue pretreatment step
- Antibody diluent, e.g. Q Diluent for IHC (Cat. No. AD-001-xxxx)
- Antigen retrieval reagent, e.g. Q Retrieval Low pH 6.0 (Cat. No. AR-001-0120) or Q Retrieval High pH 9.0 (Cat. No. AR-002-0120)
- Detection system, e.g. PolyQ Stain kits and appropriate chromogen
- Wash buffer, e.g. TBS (Cat. No. BU-006-xxxx) or TBS-Tween20 (Cat. No. BU-007-xxxx)
- Blocking reagent
- Hematoxylin
- Mounting medium
- Light microscope

Storage and handling

Store at 2 – 8 °C.

When stored correctly, the antibody is stable up to the expiration date indicated on the vial. This also applies to the shelf life after opening or after dilution of the concentrate by the end user. Do not use after expiration date.

To ensure proper reagent delivery and stability of the antibody, replace the dispenser cap after every use and immediately place the bottle cool in an upright position.

Specimen preparation

Routinely processed, FFPE tissues are suitable for use with this primary antibody. The recommended tissue fixative is 10 % neutral buffered formalin. Variable results may occur as a result of prolonged fixation or special processes such as decalcification of bone marrow preparations. Thickness of tissue sections, which should be placed on positively charged slides, should be 2 – 5 µm. Pretreatment of deparaffinized tissue with heat-induced epitope retrieval (HIER) is recommended. Slides should be stained as soon as possible, as antigenicity of cut tissue sections may diminish over time. The optimum pretreatment protocol must be determined in user's own system.

Warnings and precautions

1. Authorized and skilled personnel may only use the product.
2. There are no estimated health risks, if the product is used as directed. MSDS is available on request.
3. Product contains sodium azide as preservative. Pure sodium azide is toxic. The concentration of sodium azide in this reagent is < 0.1 % which is not classified as hazardous.
4. As with any product derived from biological sources, proper handling procedures should be used.
5. Do not use reagents after expiration date.
6. Take reasonable precautions when handling reagents. Use protective clothing and gloves.
7. All hazardous materials should be disposed according to guidelines for hazardous waste disposal. Materials of human or animal origin should be handled as biohazardous materials and disposed of with proper precautions.
8. Avoid microbial contamination of reagents as it may cause incorrect results.

Staining procedure

Primary antibody has been optimized for use in combination with PolyQ Stain detection systems. The following data are recommendations. Due to variation in tissue fixation and processing, as well as general lab instrument and environmental conditions, it may be necessary to adjust incubation times. The optimum protocol must be determined in user's own system.

Antigen retrieval: HIER; Boil tissue sections in Q Retrieval for approximately 30 min followed by cooling at room temperature (RT).

Incubation of primary antibody for 30 – 60 min at RT.

Staining protocol: Follow the procedure described in the instructions of the used detection system.

Recommendations for staining protocols:

1. Leica Bond
 - Dilution – 1:100
 - Pretreatment – HIER: ER2 30 min
 - Incubation – 60 min
 - Detection – Bond Polymer Refine

Quality control procedures

Positive tissue control

A positive tissue control must be run with every staining procedure performed for monitoring the correct performance of processed tissues and test reagents. Known positive tissue controls should not be utilized as an aid in determining a specific diagnosis of patient sample.

If the positive tissue controls fail to demonstrate appropriate positive staining, results with the test specimens must be considered invalid.

Example for positive tissue control:

- Kidney (At least weak to moderate, distinct nuclear staining should be seen in the majority of epithelial cells of the proximal and distal renal tubules, collection ducts, and the parietal epithelial cells of the Bowman's capsule)
- Fallopian tube (Weak to moderate, distinct nuclear staining in the majority of the ciliated epithelial cells and a strong nuclear staining of the intercalated secretory epithelial cells)

Negative tissue control

Negative tissue controls provide an indication of non-specific background staining. If specific staining occurs in the negative tissue control sites, results with the patient specimens must be considered invalid.

The variety of cell types present in most tissue sections offers internal negative control sites. Therefore, the same tissue used for the positive tissue control may be used as the negative tissue control.

Example for internal negative tissue control:

- Tonsil (No staining should be seen in squamous epithelial cells and in lymphocytes)

Discrepancies

If quality control results do not meet specifications, patient results are invalid. Identify and correct the problem (see section "Troubleshooting"), then repeat the entire procedure with the patient samples.

Negative control reagent

A negative control reagent is used in place of the primary antibody to evaluate non-specific staining. Host species and incubation time should be similar to primary antibody.

Interpretation of results

The immunostaining procedure causes a colored reaction product to precipitate at the antigen sites localized by the primary antibody.

Cellular localization: Nuclear.

A qualified pathologist experienced in immunohistochemistry procedures must evaluate positive and negative tissue controls before interpreting patient specimens.

Positive staining intensity should be assessed within the context of any background staining of the negative reagent control.

Note: A negative result means that the antigen in question was not detected, not that the antigen is absent in the cells or tissue assayed. A panel of antibodies may be used to verify the results. Additionally, the morphology of each tissue sample should be examined utilizing a hematoxylin and eosin stained section. A qualified pathologist must interpret the patient's morphologic findings and pertinent clinical data.

Performance characteristics

The antibody has been validated by IHC using FFPE human tissue sections of different types of healthy and neoplastic tissues.

Table 1 Testing of healthy FFPE tissue sections

Tissue	Positive/total cases
Thyroid	2/2
Kidney	10/10

Endometrium	1/1
Fallopian tube	3/3
Tonsil	0/10
Appendix	0/10

Table 2 Testing of neoplastic FFPE tissue sections

Tissue	Positive/total cases
Thyroid carcinoma	13/14
Renal cell carcinoma	9/10
Ovary carcinoma	18/20
Endometrium carcinoma	8/9

PAX8 (QR016) shows no cross reactivity with PAX5 or PAX6.

Analytical performance

The antibody passed all analytical performance tests. Analytical sensitivity has been determined by measuring concordances with known positive tissues, and analytical specificity has been determined by measuring concordances with known negative tissues with at least 90% overall concordance between the new test and the expected results each. Results are 100% each. Trueness has been verified by comparison with an independent product, and is confirmed as the results match 100%.

The method has a high level of precision - repeatability within run (carried out several times on the same day by the same analyst with the same instrument), reproducibility between run (carried out on different days by different analysts with different instruments) and reproducibility from lot to lot (results of a new reagent lot compared to those of a previously used lot) are confirmed with 100% each. Trueness and precision result in a high level of accuracy of measurement for the method.

Since IHC is a qualitative detection method that provides information on whether the corresponding antigen is present or not, the concentration-related parameters detection and quantification limits, cutoff/tolerance limit, measuring range and linearity are not applicable and are all unable to be defined for this product.

Clinical performance

For the assessment of clinical performance, authorized staining pictures provided by customers using clone QR016 have been evaluated to prove clinical evidence with the aim to achieve a robust and correct visualization of the target antigen in clinical samples with unknown expression levels, thereby contributing to a valid diagnosis.

Clinical performance parameters have been calculated with a total of 14 thyroid carcinoma, 10 renal cell carcinoma, 20 ovarian carcinoma, 9 endometrial carcinoma^[own data] as well as 10 tonsil and 10 appendix tissues^[own data]. Summarized diagnostic sensitivity, diagnostic specificity, PPV, NPV and LR- resulted in 100%, 100%, 100%, 100% and 0 (interpreted as convincing diagnostic evidence), respectively. Results of clinical performance evaluation are in accordance with scientific peer-reviewed literature as shown in section 'Summary and explanation'. The remaining tested tissues may give an overview of staining features of PAX8 (QR016).

Limitations

1. For *in vitro* diagnostic use.
2. For laboratory use only.
3. This reagent is "for professional use only" as immunohistochemistry is a multiple step process that requires specialized training in the selection of the appropriate reagents, tissues, fixation and processing, preparation of the immunohistochemistry slide, choice of detection system, and interpretation of the staining results.
4. Tissue staining is dependent on the handling, processing and storage of the tissue prior to staining.

5. Improper fixation, freezing, thawing, washing, drying, heating, sectioning, or contamination with other tissues or fluids may produce artifacts, antibody trapping, or incorrect results. Optimal performance requires adequate specimen quality as well as appropriate sample preparation.
6. Excessive or incomplete counterstaining may compromise proper interpretation of results.
7. False positive results may be seen because of non-immunological binding of proteins or substrate reaction products. They may also be caused by pseudo peroxidase activity (erythrocytes), endogenous biotin (example: liver, brain, kidney) or endogenous peroxidase activity (cytochrome C).
8. When used in blocking steps, normal sera from the same animal source as the secondary antisera may cause false negative or false positive results because of the effect of autoantibodies or natural antibodies.
9. Tissues from persons infected with hepatitis B virus and containing hepatitis B surface antigen may exhibit nonspecific staining with HRP.
10. Unexpected results may occur due to biological variability of antigen expression in neoplasms or other pathological tissues.
11. The clinical interpretation of any test results should be evaluated within the context of the patient's medical history and other diagnostic laboratory test results. Staining must be performed in a certified, licensed laboratory under the supervision of a qualified pathologist who is responsible for evaluation and assuring the adequacy of positive and negative controls. Manufacturer is not liable for incorrect results due to visual evaluation.
12. Prediluted antibodies are ready-to-use and optimized for staining. Further dilution may lead to incorrect results.
13. After successful validation users may dilute antibody concentrates according to requirements. Appropriate controls must be employed and documented.
14. The performance of the product was established using the procedures provided in this package insert only and modifications to these procedures may lead to changes in efficiency. Non-application as prescribed in this data sheet leads to loss of all liability. Any changes in product, composition, implementation, as well as use in combination with any reagents other than recommended herein is not allowed; users are responsible themselves for those changes and have to perform prior validation.
15. Application in combination with diagnostic devices, e.g. an automated staining platform, requires prior validation before staining patient specimen.
16. We do not take responsibility for any possible damage including personal injury, time or effort on economic loss caused by this product. Our warranty is limited to the price paid for the product.

Troubleshooting

1. Only intact cells should be used for interpretation of staining results, as degenerated cells show non-specific staining.
2. If no staining occurs, control application order of reagents. Follow all indications given in the instructions for use.
3. Do not allow the sections to dry out.
4. If weak staining occurs, pay attention during staining steps to freshly prepared chromogen, incubation times and temperatures, as well as accurate draining off of reagents.
5. Avoid surplus background staining by optimal removal of paraffin, washing of slides and dilution of primary antibody. If excessive background staining occurs, high levels of endogenous biotin may be present (unless a biotin-free detection system is being used). A biotin blocking step should be included.

6. Sodium azide inactivates HRP, which may lead to false results. Wash sections in sodium azide free buffer.
7. Contact quartett customer service in case of any uncertainties.

Literature

[1] Tacha D, Zhou D, Cheng L (2011). Appl Immunohistochem Mol Morphol. 19(4):293-9.

[2] Laury AR, Perets R, Piao H et al. (2011). Am J Surg Pathol.35(6):816-26.

Distributor

quartett Biotechnologie GmbH
 Am Mühlenberg 4, 14476 Potsdam, Germany
 Tel: +49 (0)30 765 925-0 • Fax: +49 (0)30 765 925-55
 service@quartett.com • www.quartett.com

Manufacturer

biocyc Biotechnologie GmbH & Co. KG
 Am Mühlenberg 4, 14476 Potsdam, Germany
 cert. by TÜV Rheinland LGA Products GmbH
 ISO 13485 & ISO 9001
 Tel: +49 (0)331 967 826-00

In the event that the user experiences any technical or performance-related issues with the product, please consult the manufacturer or a competent authority.


Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the member state in which the user and/or the patient is established.

Summary of safety and performance (SSP) can be found in EUDAMED when the related module is available.

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