AH/AHS-B Range Automatic Autoclaves with fractioned vacuum

AH-21 B

AHS-50 B


AH/AHS-B Range

B Class: Suitable to sterilize: unwrapped and wrapped solids, porous objects, and hollow ware.

Operating Principal

The purpose of the sterilization with pre-vacuum is the efficient steam penetration throughout the entire load. This performance is achieved with a vacuum pump, which removes air by fractioned pre-vacuum.

A built-in cooling system by radiator coil enables a better performance of the vacuum pump, removing air pockets from all load types.

Immediate production of steam through the built-in steam generator allows optimization of cycle times.

Touch Screen Display – Brand new Touch Screen color display for an easy programming and parameters selection for sterilization process.

AH/AHS-B Range

Main screen


Sterilization screen

Programs – Easy to operate with 12 pre-defined programs

Nº Program	Program name	Nº Pre-vacuum	Temperature (°C)	Sterilization time (min.)	Drying time (min.)
P1	Helix-Bowie&Dick Test	3	134	4'	4'
P2	Vacuum Test	1	-	-	-
P3	Porous-134	3	134	4'	10'
P4	Prion-134	3	134	18'	10'
P5	Porous-121	3	121	20'	10'
P6	Hollowware-134°C	3	134	4'	5'
P7	Hollowware -121°C	3	121	20'	5'
P8	Wrapped-134°C	1	134	4'	5'
P9	Wrapped -121°C	1	121	20'	5'
P10	Solids-134°C	1	134	4'	5'
P11	Solids -121°C	1	121	20'	5'
P12	Flash-134°C	1	134	3'	1'

Up to 50 programs user-defined with the following parameters setting:

- Sterilization temperature
- Sterilization time
- Drying time
- Number of pre-vacuum

Parameter selection screen

AH/AHS-B Range

Water Filling – Pure water for sterilization can be added manually into the water pan. Also, automatic filling is available by connecting the device to pure mains water.

Sterilization Process

Heating Stage – Fast saturated steam production by means of a powerful built-in steam generator.

Pre-vacuum Stage – Vacuum pump with built-in cooling system by radiator coil provides a better vacuum performance.

Sterilization Control – A PT-100 temperature probe controls the chamber temperature, with the possibility of installing another flexible PT-100 temperature probe for liquid programs.

Drying Stage – Powerful heating jacket together with fractionated pre-vacuum provides a fast and intense drying.

Security

All safety standards are covered with safety valve, safety thermostat and safety pressure switch, door opening prevention system in case of overpressure, open door sensor and water level detector. EU Directive and Standards:

2014/35/UE *Low Voltage*

2014/30/UE *Electromagnetic Compatibility*

2014/68/UE *Pressure Equipment*

EN-61010-1

EN-61010-2-040

EN-61326

AD 2000 Merkblatt

Communications and Documentation

USB Port – For electronic storage of process data to a USB flash drive.


Ethernet Port – For live monitoring with direct data storage. It is needed the customized software SW8000 (optional).

Integrated Printer (optional) – For quality printouts after each cycle. Intervals can be adjustable. Printer has to be requested before ordering.

AH/AHS-B Range


Components supplied with the autoclave AH-21 B

1. Tray support for four trays.
2. Three stainless steel perforated trays. Ref. BAH-21
3. Auxiliary tray.
4. A holding clamp.
5. 1.5 m flexible tube to drain condensation reservoir.
6. Bag support. Ref. BAP-21

Components supplied with the autoclaves AHS-50/75 B

1. Support for five trays.
2. Two perforated trays made of stainless steel. Ref. BAH-B
3. Auxiliary tray.
4. 1.5 m flexible tube to drain condensation reservoir.

AH/AHS-B Range

Technical Data

Referencia Reference Référence	Volumen total / útil Total volume / usable Volume total / utile l	Dimensiones útiles cámara Usable chamber dimensions Dimensions utile chambre H max. x D mm	EXTERIOR L x D x H mm	Temp. °C	Presión máxima Max. pressure Pression max. bar	Potencia Power Puissance W	Peso bruto Gross weight Poids brut Kg	Tensión Voltage Tension V	Frecuencia Frequency Fréquence Hz	Depósito / Container / Réservoir l	
										AGUA LIMPIA CLEAN WATER EAU PURE	AGUA USADA WASTE WATER EAU UTILISÉE
AH-21-B	22/21	210 x 430	560 x 660 x 425	105-134	2,1	2300	85	230	50/60	6	2,6
AHS-50-B	55/50	360 x 400	805 x 805 x 650	105-134	2,1	3600	125	230	50/60	10	6
AHS-75-B	79/75	360 x 600	805 x 1005 x 650	105-134	2,1	3600	140	230	50/60	10	6

Consultar por otras tensiones / Ask for other voltages / Demander d'autres tensions

AHS-50/75 models with closing wheel

Electrical connection for equipment of 21, 50 and 75 L: 230V (1P + N), 50/60 Hz.

Audible and visible end-of-cycle sterilization indicator

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FLEXIBLE "HEART" PROBE PT-100

- Recommended for sterilizing liquids in flasks of 1 litre or bigger.
- The control is done from the microprocessor.
- This accessory should be installed in factory.

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THERMAL PAPER PRINTER

- It prints program number, cycle number, temperature, date and hour of the run.
- Error messages are also printed.
- To be installed in our facilities.

References:
IT: LCD Screen
IT/TS: Touch Screen

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SET AUTOMATIC FILLING


- To be installed in our facilities.

References:
KLL: for B and DRY from 28 up to 150 L
KLL-21: for AH-21B and AH-21L
KLL-MP: for MP

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* Serie AE-B, AES-RFG, AE-MP

PC CONNECTION SW8000

- Communication software between the equipment and the PC.
- Connection for Ethernet.
- USB Adapter.
- Compatible with touch screen circuits.

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* Todos los modelos / All models / Tous les modèles

STERILIZATION CONTROL TAPE

- Chemical indicator for steam. The change of color indicates that the materials have been processed, without being a guarantee of proper sterilization. Class 1.
- To verify that sterilization is effective, additional methods are needed as biological indicators EN ISO 11138.
- Tape roll of 50 m x 19 mm

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* Serie AHS, AH-B, AE-B

BOWIE & DICK TEST PACK

- Class B indicator.
- Check the correct steam penetration in the porous bodies.
- Indicator printed with non-toxic inks and laminated.
- Box of 30 tests.

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* Serie AH-B

HELIX PCD TEST

- Class B indicator.
- Check the correct steam penetration in the close lumen bodies.

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