Free T3

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 Free T3

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of free T3 in serum. Measurements of free T3 are used in the diagnosis and treatment of thyroid disease.

Catalog Numbers:
L2KF32 (200 tests)
L2KF36 (600 tests)

Test Code: **FT3** Color: **Aqua**

Summary and Explanation

Under normal physiological conditions triiodothyronine (T3) represents approximately 5 percent of the thyroid hormone in plasma. Although present in lower concentration than thyroxine (T4), T3 has greater metabolic activity, faster turnover and a larger volume of distribution. It is produced largely through extrathyroidal conversion of T4. Like T4, it circulates almost entirely bound to the carrier proteins TBG, pre-albumin and albumin. Free T3 constitutes only about 0.25% of the total T3 in circulation.

Measurement of total T3 by immunoassay has a number of well-established uses. In the presence of elevated free or total T4, total T3 measurements help confirm a hyperthyroidism diagnosis. Abnormal elevation of total T3 may also occur when the total T4 concentration is normal — a condition known as "T3 toxicosis".

For the most part, free T3 levels correlate closely with total T3 levels. Total T3, however, depends not only on thyroid status and the peripheral conversion of T4 to T3, but also on the concentration of thyroid hormone-binding proteins. Free T3, on the other hand, is largely unaffected by variations in these carrier proteins. Thus, the TBG elevations typical of pregnancy, oral contraceptive use, and estrogen therapy effects an increase in the total T3 level while leaving the free T3 concentration basically unchanged.

The free T3 concentration typically reflects a patient's actual thyroid status more reliably than the total T3 concentration.

Principle of the Procedure

Competitive, analog-based immunoassay.^{25,26}

Incubation Cycles: 2 × 30 minutes

Specimen Collection

Because EDTA would affect results, it should not be used as an anticoagulant.

Lipemic or grossly contaminated samples may give erroneous results. The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Free T3 has not been tested with all possible variations of tube types.

Volume Required: 100 µL serum

Storage: 2 days at 2–8°C, or 2 months at –20°C.

Warnings and Precautions

For *in vitro* diagnostic use.

CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.²⁷⁻²⁹

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. Contains: sodium azide; Free T3 Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Free T3 Bead Pack (L2F312)

With barcode. 200 beads, coated with monoclonal murine anti-T3. Stable at 2–8°C until expiration date.

L2KF32: 1 pack

L2KF36: 3 packs

Free T3 Reagent Wedge (L2F3A2)

With barcode, containing three liquid reagents, 10.5 mL, 10.5 mL and 11.5 mL respectively, ready to use. The first and second reagents consist of ligand-labeled T3 analog in buffer, with preservative. The third reagent consists of alkaline phosphatase (bovine calf intestine) conjugated to anti-ligand in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KF32: 1 wedge

L2KF36: 3 wedges

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

Free T3 Adjustors (LF3L, LF3H)

2 vials (Low and High), of lyophilized free T3 in processed human serum, with preservative. Reconstitute each vial with 4.0 mL distilled or deionized water. Stable at 2–8°C for 30 days after reconstitution, or for 6 months (aliquotted) at –20°C.

L2KF32: 1 set

L2KF36: 2 sets

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

Also Required

Distilled or deionized water; test tubes; controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of free T3.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Expected Values

Adult: A reference range study for IMMULITE 2000 Free T3 was performed on 212 apparently healthy adults yielding a central 95% range of 1.8–4.2 pg/mL for euthyroid adults.

Pediatric: Reference intervals for the pediatric population (children and adolescents) were established for the IMMULITE Free T3 assay in accordance with CLSI guideline EP28-A3C.³⁰

For analysis of data, the population was divided into three age subgroups:

- Infants: subjects aged 1–23 months
- Children: subjects aged 2–12 years
- Adolescents: subjects aged 13–20 years

A non-parametric approach was used to establish the reference intervals for children and adolescents where the 2.5 and 97.5 percentiles of the distribution of values were calculated. For the infant population, a robust measure of location and spread, as developed by Horn and

Pesce, was used to estimate the 2.5 and the 97.5 percentile reference intervals, accommodating the smaller sample size.^{30–32}

The reference intervals detailed by age group and number of samples are presented in the Reference Intervals table.

IMMULITE 2000 Free T3 Pediatric Reference Intervals

Age Group	<i>n</i>	Conventional (pg/mL)	SI (pmol/L)
Infants (1–23 Months)	81	3.6–7.5	5.5–11.5
Children (2–12 Years)	195	3.7–6.6	5.7–10.1
Adolescents (13–20 Years)	148	3.1–5.9	4.8–9.1

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

The free T3 assay serves primarily as a confirmatory test for hyperthyroidism and is of very limited value in the context of hypothyroidism. There may be some overlap between the euthyroid and hyperthyroid ranges.

The interpretation of free T3 results is complicated by a variety of drugs which can affect the binding of T3 to the thyroid hormone carrier proteins.

In severe nonthyroidal illness, the assessment of thyroid status becomes especially difficult. Since some patients in this category may suffer from concomitant primary hypothyroidism or from a compensatory secondary hypothyroidism, TSH immunoassay has been recommended as a confirmatory test in this context.

In rare conditions associated with extreme variations in the albumin-binding capacity — such as familial dysalbuminemic hyperthyroxinemia (FDH) — direct free thyroid hormone assays may yield misleading results.

Circulating autoantibodies to T3, and hormone-binding inhibitors, may interfere with the assay.^{2,3,12}

Heparin has been reported to have both *in vivo* and *in vitro* effects on free thyroid hormones.⁸ Hence, samples should not be collected during or soon after the administration of this anticoagulant.

Since dilution shifts the equilibrium between free and protein-bound T3, the assay cannot be expected to maintain linearity under dilution. Accordingly, no attempt should be made to dilute samples with high free T3 results. In particular, any sample with a result greater than the upper limit of the assay's calibration range should simply be reported as such.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in pg/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Conversion Factor:

pg/mL $\times$ 1.536 $\rightarrow$ pmol/L

Calibration Range: 1.0–40 pg/mL
(1.5–61 pmol/L)

The assay is traceable to an internal standard manufactured using qualified materials and measurement procedures.

Analytical Sensitivity: 1.0 pg/mL
(1.5 pmol/L)

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

Specificity: The antibody is highly specific for free T3. (See "Specificity" table.)

Effect of TBG: The assay's zero calibrator was spiked with T3-stripped TBG resin (20, 40 and 80 μ g/mL) and assayed. The results showed no interference. (See "Effect of TBG" table.)

Effect of Albumin: Samples spiked with charcoal-absorbed human serum albumin (10, 20 and 50 mg/mL) were assayed. (See "Effect of Albumin" table.)

Effect of Nonesterified Fatty Acids

(NEFA): Samples spiked with oleic acid (2.5 and 5.0 mmol/L) were assayed. Results show no significant effect at 2.5 mmol/L, but decreased values at the higher level. (See "Effect of Nonesterified Fatty Acids" table. Increased levels of NEFA may be encountered in late pregnancy and in heparin therapy.)

Other Compounds: Results for samples assayed after spiking with phenylbutazone (10, 100 and 1000 μ g/dL), or phenytoin (5, 10, 25, and 50 μ g/mL), or salicylate (10, 25, 50 and 100 mg/dL) showed no effect from these compounds. (See "Effect of Phenylbutazone", "Effect of Phenytoin" and "Effect of Salicylate" tables.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Biotin: Specimens that contain biotin at a concentration of 1500 ng/mL demonstrate a less than or equal to 10% change in results. Biotin concentrations greater than this may lead to incorrect results for patient samples.

Hemolysis: Presence of hemoglobin in concentrations up to 384 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Highly lipemic samples may interfere with the assay. (See "Effect of Lipemia" table.)

Method Comparison: The assay was compared to IMMULITE Free T3 on 105 samples. (Concentration range: approximately 1 to 14 pg/mL.)
By linear regression:

(IML 2000) = 1.05 (IML) + 0.41 pg/mL
 $r = 0.977$

Means:

3.5 pg/mL (IMMULITE 2000)

3.0 pg/mL (IMMULITE)

References

- 1) Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982.
- 2) Beck-Peccoz P, Romelli PB, et al. Free T4 and free T3 measurement in patients with anti-iodothyronine autoantibodies. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 231-8.
- 3) Bhagat CI, et al. Interference in radioimmunoassay of total serum thyroxine and free thyroxine due to thyroxine-binding autoantibodies. Clin Chem 1983;29:1324-5.
- 4) Borst GC, Eil C, Burman KD. Euthyroid hyperthyroxinemia; a review. Ann Intern Med 1983;98:366-78.
- 5) Burger AG, Lambert MJ, Cullen M. Pharmacological influence on free thyroid hormones. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 303-14.
- 6) Caroyan P, et al. Assessment and clinical significance of free thyroid hormone radioimmunoassays. In: Ekins RP, et al, editors. Free thyroid hormones. Amsterdam: Excerpta Medica, 1979: 181-93.
- 7) Carayon P, et al. Free thyroid hormones in euthyroidism: evaluation of normal ranges in various physiological conditions. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 179-85.
- 8) Diagnostic Products Corporation. Coat-A-Count Free T4, package insert (PITKF4).
- 9) Ekins RP. Methods for the measurement of free thyroid hormones. In: Ekins RP, editor. Free thyroid hormones. Amsterdam: Excerpta Medica, 1979: 72-106.
- 10) Ekins RP. The radioimmunoassay of free thyroid hormones in blood. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 73-90.
- 11) Ellis SM, Ekins RP. The radioimmunoassay of serum free triiodothyronine and thyroxine. In: Pasternak CA, editor. Radioimmunoassay in clinical biochemistry. London: Heyden, 1975: 187-94.
- 12) John R. Autoantibodies to thyroxine and interference with free thyroxine assay. Clin Chem 1983;29:581-2.
- 13) Larsen PR. The use of serum T3 measurements by radioimmunoassay in the diagnosis of thyroid disease. In: Abraham GE, editor. Radioassay systems in clinical endocrinology. New York: Marcel Dekker, 1981: 117-29.
- 14) Parlow ME, Oddie TH, Fisher DA. Evaluation of serum triiodothyronine and adjusted triiodothyronine (free triiodothyronine index) in pregnancy. Clin Chem 1977;23:490-2.
- 15) Lin T, Nankin HR. Direct measurement of free thyroxine in patients with T3-toxicosis. Horm Metab Res 1980;12:611-4.
- 16) Mardell R, Gerson M. A method of assessing serum triiodothyronine concentrations that is independent of subject's age and variations in concentrations of binding proteins in serum. Clin Chem 1978;24:1792-6.
- 17) Pennisi F, Romelli PB, Vancheri L. Measurement of free thyroid hormones in serum by column adsorption chromatography and radioimmunoassay. In: Albertini A, Ekins RP, editors. Free hormones in blood. Amsterdam: Elsevier Biomedical Press, 1982: 107-12.
- 18) Roux F, et al. Clinical merits in the radioimmunoassay of free serum tetraiodothyronine (FT4) and free serum triiodothyronine (FT3). Int J Nuclear Med Biol 1980;7:386-90.
- 19) Salmon D, et al. Chemical hyperthyroidism: serum triiodothyronine levels in clinically euthyroid individuals treated with levothyroxine. Arch Intern Med 1982;142:571-3.
- 20) Sawin CT, Chopra D, et al. The free triiodothyronine (T3) index. Ann Intern Med 1978;88:474-7.
- 21) Siegel L, McDonald LJ, Robin NI. Estimation of free triiodothyronine in serum: a new method and its clinical relevance. Clin Chem 1978;24:1891-4.
- 22) Smals AGH, Ross AH, Kloppenborg PWC. Dichotomy between serum free triiodothyronine and free thyroxine concentrations in familial thyroxine-binding globulin deficiency. J Clin Endocrinol Metab 1981;53:917-22.
- 23) Walfish PG. Thyroid physiology and pathology. In: Collu R, Ducharme JR, Guyda H, editors. Pediatric endocrinology. New York: Raven Press, 1981: 357-431.
- 24) Weeke J, et al. A longitudinal study of serum TSH, and total and free iodothyronines during normal pregnancy. Acta Endocrinol 1982;101:531-7.
- 25) Witherspoon LR, El Shami AS, et al. Chemically blocked analog assays for free thyronines. Clin Chem 1988;34:9-16 and 17-23.
- 26) Diagnostic Products Corporation. Albumin-related artifacts in free thyroid hormone assays. Los Angeles, 1987 (Document number: ZJ026).
- 27) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8.
- 28) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
- 29) NCCLS Document M29-A3.
- 30) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.
- 31) Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Guideline EP28-A3C.
- 32) Horn PS, Pesce AJ, Reference Intervals. A User's Guide, Washington, DC: AACC Press; 2005.
- 33) Reed AH, Henry RJ, Mason WB. Influence of statistical method used on the resulting estimate of normal range. 1971;17:275-284.

Technical Assistance

In the United States, contact
Siemens Healthcare Diagnostics
Technical Services department.
Tel: 877.229.3711. Outside the United
States, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare
Diagnostics Products Ltd. is certified to ISO
13485.

Tables and Graphs


Precision (pg/mL)

Mean ³	Within-Run ¹			Total ²	
	SD ⁴	CV ⁵	SD	SD	CV
1	2.5	0.21	8.4%	0.25	10%
2	3.2	0.29	9.1%	0.30	9.4%
3	3.8	0.34	9.0%	0.38	10%
4	4.2	0.24	5.7%	0.34	8.1%
5	5.7	0.31	5.4%	0.40	7.0%
6	6.5	0.38	5.9%	0.48	7.4%
7	13	0.56	4.3%	0.70	5.4%

Effect of Albumin (pg/mL)

Unspiked ²	Albumin (mg/mL Added) ¹			
	10	20	50	
1	2.0	1.8	1.9	2.4
2	2.0	2.3	2.0	1.7
3	2.9	2.9	2.9	3.0
4	3.0	3.0	3.5	3.1
5	3.0	3.1	3.4	3.0
6	3.8	3.9	4.1	3.8

Method Comparison

$$(IML\ 2000) = 1.05 (IML) + 0.41\text{ pg/mL}$$
$$r = 0.977$$

Effect of TBG (pg/mL)

TBG ($\mu\text{g/mL}$ Added) ¹	Binding (%B/B ₀) ²
0	100%
20	92%
40	95%
80	90%

Effect of Nonesterified Fatty Acids (pg/mL)

Unspiked ²	Oleic Acid (mmol/L Added) ¹		
	2.5	5.0	
1	2.1	2.6	1.3
2	2.4	2.9	1.7
3	2.8	2.9	1.5
4	2.9	2.7	1.6
5	3.0	3.2	1.9
6	3.5	3.8	2.4

Effect of Phenylbutazone (pg/mL)

Unspiked ²	Phenylbutazone ($\mu\text{g/dL}$ Added) ¹			
	10	100	1000	
1	2.2	2.0	2.3	2.5
2	2.5	2.4	2.7	2.8
3	2.6	2.6	2.8	2.6
4	3.1	3.3	3.2	3.2
5	3.3	3.3	3.3	3.2
6	3.6	3.4	3.2	3.7

Effect of Salicylate (pg/mL)

		Salicylate (mg/dL Added) ¹			
	Unspiked ²	10	25	50	100
1	2.2	2.9	2.7	2.9	2.8
2	2.7	2.9	3.0	2.6	2.3
3	2.8	2.6	2.4	2.5	2.6
4	3.6	3.2	4.1	3.8	3.0
5	4.0	3.7	3.9	3.7	3.5
6	4.0	4.5	4.3	4.1	3.7
7	4.2	4.5	4.2	4.0	3.3
8	4.4	3.7	4.0	3.3	3.5
9	6.9	8.0	7.7	7.6	6.7
10	9.7	9.4	9.3	8.7	8.8
11	12.4	12.2	12.0	11.0	10.8
12	16	15	15	15	14

Effect of Phenytoin (pg/mL)

		Phenytoin (µg/mL Added) ¹			
	Unspiked ¹	5	10	25	50
1	2.5	2.5	2.6	2.5	2.8
2	2.9	2.6	2.5	2.6	2.7
3	3.0	2.9	2.5	2.7	2.7
4	3.0	2.9	2.9	3.2	3.1
5	6.4	6.8	6.9	7.2	7.1
6	7.7	7.6	7.5	7.0	7.6

Effect of Lipemia

	Triglycerides Added mg/dL	Observed ¹	Expected ²	%O/E ³
1	—	2.0		
	250	2.1	2.0	105%
	500	2.1	1.9	111%
	1000	1.8	1.9	95%
	2000	1.7	1.8	94%
	3000	1.6	1.7	94%
2	—	2.1		
	250	1.8	2.1	86%
	500	2.1	2.1	100%
	1000	1.9	2.0	95%
	2000	1.9	1.9	100%
	3000	2.1	1.8	117%
3	—	2.6		
	250	2.4	2.6	92%
	500	2.3	2.6	88%
	1000	2.4	2.5	96%
	2000	1.7	2.4	71%
	3000	1.4	2.2	64%
4	—	2.7		
	250	2.7	2.7	100%
	500	2.3	2.7	85%
	1000	2.5	2.6	96%
	2000	2.2	2.5	88%
	3000	2.3	2.3	100%
5	—	2.9		
	250	2.8	2.9	97%
	500	2.6	2.8	93%
	1000	2.3	2.8	82%
	2000	2.3	2.6	88%
	3000	1.9	2.5	76%
6	—	3.3		
	250	2.8	3.3	85%
	500	3.3	3.2	103%
	1000	2.7	3.1	87%
	2000	2.6	3.0	87%
	3000	2.4	2.8	86%

Specificity

Compound ¹	pg/mL Added ²	% Cross- reactivity ³
L-Thyroxine (T4)	50,000	ND
D-Thyroxine	50,000	ND
Triiodothyroacetic acid	5000	0.073%
Tetraiodothyroacetic acid	200,000	ND
Diiodo-L-tyrosine	10,000,000	ND
Methimazole	10,000,000	0.00002%
5,5'-Diphenylhydantoin	10,000,000	ND
Phenylbutazone	10,000,000	ND
6- <i>n</i> -Propyl-2-thiouracil	10,000,000	ND
Triiodo-L-tyrosine	10,000,000	0.00002%

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IMMULITE and Coat-A-Count are trademarks of Siemens Healthcare Diagnostics.
All rights reserved.